



US 20050042606A9

(19) **United States**
(12) **Patent Application Publication**
Bergeron et al.

(10) **Pub. No.: US 2005/0042606 A9**
(48) **Pub. Date: Feb. 24, 2005**
CORRECTED PUBLICATION

(54) **SPECIFIC AND UNIVERSAL PROBES AND AMPLIFICATION PRIMERS TO RAPIDLY DETECT AND IDENTIFY COMMON BACTERIAL PATHOGENS AND ANTIBIOTIC RESISTANCE GENES FROM CLINICAL SPECIMENS FOR ROUTINE DIAGNOSIS IN MICROBIOLOGY LABORATORIES**

(76) Inventors: **Michel G. Bergeron**, Sillery (CA); **Marc Ouellette**, Quebec (CA); **Paul H. Roy**, Loretteville (CA)

Correspondence Address:
QUARLES & BRADY LLP
411 E. WISCONSIN AVENUE
SUITE 2040
MILWAUKEE, WI 53202-4497 (US)

(21) Appl. No.: **10/121,120**
(22) Filed: **Apr. 11, 2002**

Prior Publication Data

(15) Correction of US 2003/0180733 A1 Sep. 25, 2003
See Related U.S. Application Data.
(65) US 2003/0180733 A1 Sep. 25, 2003

Related U.S. Application Data

(63) Continuation of application No. 09/452,599, filed on Dec. 1, 1999, now abandoned.
Continuation of application No. 08/526,840, filed on Sep. 11, 1995, now Pat. No. 6,001,564, which is a continuation-in-part of application No. 08/304,732, filed on Sep. 12, 1994, now abandoned.

Publication Classification

(51) **Int. Cl.⁷** **C12Q 1/68**; C12P 19/34
(52) **U.S. Cl.** **435/6**; 435/91.2

(57) **ABSTRACT**

The present invention relates to DNA-based methods for universal bacterial detection, for specific detection of the common bacterial pathogens *Escherichia coli*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa*, *Proteus mirabilis*, *Streptococcus pneumoniae*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Enterococcus faecalis*, *Staphylococcus saprophyticus*, *Streptococcus pyogenes*, *Haemophilus influenzae* and *Moraxella catarrhalis* as well as for specific detection of commonly encountered and clinically relevant bacterial antibiotic resistance genes directly from clinical specimens or, alternatively, from a bacterial colony. The above bacterial species can account for as much as 80% of bacterial pathogens isolated in routine microbiology laboratories.

The core of this invention consists primarily of the DNA sequences from all species-specific genomic DNA fragments selected by hybridization from genomic libraries or, alternatively, selected from data banks as well as any oligonucleotide sequences derived from these sequences which can be used as probes or amplification primers for PCR or any other nucleic acid amplification methods. This invention also includes DNA sequences from the selected clinically relevant antibiotic resistance genes.

With these methods, bacteria can be detected (universal primers and/or probes) and identified (species-specific primers and/or probes) directly from the clinical specimens or from an isolated bacterial colony. Bacteria are further evaluated for their putative susceptibility to antibiotics by resistance gene detection (antibiotic resistance gene specific primers and/or probes). Diagnostic kits for the detection of the presence, for the bacterial identification of the above-mentioned bacterial species and for the detection of antibiotic resistance genes are also claimed. These kits for the rapid (one hour or less) and accurate diagnosis of bacterial infections and antibiotic resistance will gradually replace conventional methods currently used in clinical microbiology laboratories for routine diagnosis. They should provide tools to clinicians to help prescribe promptly optimal treatments when necessary. Consequently, these tests should contribute to saving human lives, rationalizing treatment, reducing the development of antibiotic resistance and avoid unnecessary hospitalizations.

SPECIFIC AND UNIVERSAL PROBES AND AMPLIFICATION PRIMERS TO RAPIDLY DETECT AND IDENTIFY COMMON BACTERIAL PATHOGENS AND ANTIBIOTIC RESISTANCE GENES FROM CLINICAL SPECIMENS FOR ROUTINE DIAGNOSIS IN MICROBIOLOGY LABORATORIES

BACKGROUND OF THE INVENTION

[0001] Classical Identification of Bacteria

[0002] Bacteria are classically identified by their ability to utilize different substrates as a source of carbon and nitrogen through the use of biochemical tests such as the API20E™ system. Susceptibility testing of Gram negative bacilli has progressed to microdilution tests. Although the API and the microdilution systems are cost-effective, at least two days are required to obtain preliminary results due to the necessity of two successive overnight incubations to isolate and identify the bacteria from the specimen. Some faster detection methods with sophisticated and expensive apparatus have been developed. For example, the fastest identification system, the autoSCAN-Walk-Away system™ identifies both Gram negative and Gram positive from isolated bacterial colonies in 2 hours and susceptibility patterns to antibiotics in only 7 hours. However, this system has an unacceptable margin of error, especially with bacterial species other than *Enterobacteriaceae* (York et al., 1992. J. Clin. Microbiol. 30:2903-2910). Nevertheless, even this fastest method requires primary isolation of the bacteria as a pure culture, a process which takes at least 18 hours if there is a pure culture or 2 to 3 days if there is a mixed culture.

[0003] Urine Specimens

[0004] A large proportion (40-50%) of specimens received in routine diagnostic microbiology laboratories for bacterial identification are urine specimens (Pezzlo, 1988, Clin. Microbiol. Rev. 1:268-280). Urinary tract infections (UTI) are extremely common and affect up to 20% of women and account for extensive morbidity and increased mortality among hospitalized patients (Johnson and Stamm, 1989; Ann. Intern. Med. 111:906-917). UTI are usually of bacterial etiology and require antimicrobial therapy. The Gram negative bacillus *Escherichia coli* is by far the most prevalent urinary pathogen and accounts for 50 to 60% of UTI (Pezzlo, 1988, op. cit.). The prevalence for bacterial pathogens isolated from urine specimens observed recently at the "Centre Hospitalier de l'Université Laval (CHUL)" is given in Tables 1 and 2.

[0005] Conventional pathogen identification in urine specimens. The search for pathogens in urine specimens is so preponderant in the routine microbiology laboratory that a myriad of tests have been developed. The gold standard is still the classical semi-quantitative plate culture method in which a calibrated loop of urine is streaked on plates and incubated for 18-24 hours. Colonies are then counted to determine the total number of colony forming units (CFU) per liter of urine. A bacterial UTI is normally associated with a bacterial count of $\geq 10^7$ CFU/L in urine. However, infections with less than 10^7 CFU/L in urine are possible, particularly in patients with a high incidence of diseases or those catheterized (Stark and Maki, 1984, N. Engl. J. Med. 311:560-564). Importantly, close to 80% of urine specimens tested are considered negative ($<10^7$ CFU/L; Table 3).

[0006] Accurate and rapid urine screening methods for bacterial pathogens would allow a faster identification of negative results and a more efficient clinical investigation of the patient. Several rapid identification methods (Uriscreen™, UTIscreen™, Flash Track™ DNA probes and others) were recently compared to slower standard biochemical methods which are based on culture of the bacterial pathogens. Although much faster, these rapid tests showed low sensitivities and specificities as well as a high number of false negative and false positive results (Koenig et al., 1992. J. Clin. Microbiol. 30:342-345; Pezzlo et al., 1992. J. Clin. Microbiol. 30:640-684).

[0007] Urine specimens found positive by culture are further characterized using standard biochemical tests to identify the bacterial pathogen and are also tested for susceptibility to antibiotics.

[0008] Any Clinical Specimens

[0009] As with urine specimen which was used here as an example, our probes and amplification primers are also applicable to any other clinical specimens. The DNA-based tests proposed in this invention are superior to standard methods currently used for routine diagnosis in terms of rapidity and accuracy. While a high percentage of urine specimens are negative, in many other clinical specimens more than 95% of cultures are negative (Table 4). These data further support the use of universal probes to screen out the negative clinical specimens. Clinical specimens from organisms other than humans (e.g. other primates, mammals, farm animals or live stocks) may also be used.

[0010] Towards the Development of Rapid DNA-Based Diagnostic Tests

[0011] A rapid diagnostic test should have a significant impact on the management of infections. For the identification of pathogens and antibiotic resistance genes in clinical samples, DNA probe and DNA amplification technologies offer several advantages over conventional methods. There is no need for subculturing, hence the organism can be detected directly in clinical samples thereby reducing the costs and time associated with isolation of pathogens. DNA-based technologies have proven to be extremely useful for specific applications in the clinical microbiology laboratory. For example, kits for the detection of fastidious organisms based on the use of hybridization probes or DNA amplification for the direct detection of pathogens in clinical specimens are commercially available (Persing et al, 1993. Diagnostic Molecular Microbiology: Principles and Applications, American Society for Microbiology, Washington, D.C.).

[0012] The present invention is an advantageous alternative to the conventional culture identification methods used in hospital clinical microbiology laboratories and in private clinics for routine diagnosis. Besides being much faster, DNA-based diagnostic tests are more accurate than standard biochemical tests presently used for diagnosis because the bacterial genotype (e.g. DNA level) is more stable than the bacterial phenotype (e.g. biochemical properties). The originality of this invention is that genomic DNA fragments (size of at least 100 base pairs) specific for 12 species of commonly encountered bacterial pathogens were selected from genomic libraries or from data banks. Amplification primers or oligonucleotide probes (both less than 100 nucleotides in

length) which are both derived from the sequence of species-specific DNA fragments identified by hybridization from genomic libraries or from selected data bank sequences are used as a basis to develop diagnostic tests. Oligonucleotide primers and probes for the detection of commonly encountered and clinically important bacterial resistance genes are also included. For example, Annexes I and II present a list of suitable oligonucleotide probes and PCR primers which were all derived from the species-specific DNA fragments selected from genomic libraries or from data bank sequences. It is clear to the individual skilled in the art that oligonucleotide sequences appropriate for the specific detection of the above bacterial species other than those listed in Annexes 1 and 2 may be derived from the species-specific fragments or from the selected data bank sequences. For example, the oligonucleotides may be shorter or longer than the ones we have chosen and may be selected anywhere else in the identified species-specific sequences or selected data bank sequences. Alternatively, the oligonucleotides may be designed for use in amplification methods other than PCR. Consequently, the core of this invention is the identification of species-specific genomic DNA fragments from bacterial genomic DNA libraries and the selection of genomic DNA fragments from data bank sequences which are used as a source of species-specific and ubiquitous oligonucleotides. Although the selection of oligonucleotides suitable for diagnostic purposes from the sequence of the species-specific fragments or from the selected data bank sequences requires much effort it is quite possible for the individual skilled in the art to derive from our fragments or selected data bank sequences suitable oligonucleotides which are different from the ones we have selected and tested as examples (Annexes I and II).

[0013] Others have developed DNA-based tests for the detection and identification of some of the bacterial pathogens for which we have identified species-specific sequences (PCT patent application Serial No. WO 93/03186). However, their strategy was based on the amplification of the highly conserved 16S rRNA gene followed by hybridization with internal species-specific oligonucleotides. The strategy from this invention is much simpler and more rapid because it allows the direct amplification of species-specific targets using oligonucleotides derived from the species-specific bacterial genomic DNA fragments.

[0014] Since a high percentage of clinical specimens are negative, oligonucleotide primers and probes were selected from the highly conserved 16S or 23S rRNA genes to detect all bacterial pathogens possibly encountered in clinical specimens in order to determine whether a clinical specimen is infected or not. This strategy allows rapid screening out of the numerous negative clinical specimens submitted for bacteriological testing.

[0015] We are also developing other DNA-based tests, to be performed simultaneously with bacterial identification, to determine rapidly the putative bacterial susceptibility to antibiotics by targeting commonly encountered and clinically relevant bacterial resistance genes. Although the sequences from the selected antibiotic resistance genes are available and have been used to develop DNA-based tests for their detection (Ehrlich and Greenberg, 1994. PCR-based Diagnostics in Infectious Diseases, Blackwell Scientific Publications, Boston, Mass.; Persing et al, 1993. Diagnostic Molecular Microbiology: Principles and Applications,

American Society for Microbiology, Washington, D.C.), our approach is innovative as it represents major improvements over current "gold standard" diagnostic methods based on culture of the bacteria because it allows the rapid identification of the presence of a specific bacterial pathogen and evaluation of its susceptibility to antibiotics directly from the clinical specimens within one hour.

[0016] We believe that the rapid and simple diagnostic tests not based on cultivation of the bacteria that we are developing will gradually replace the slow conventional bacterial identification methods presently used in hospital clinical microbiology laboratories and in private clinics. In our opinion, these rapid DNA-based diagnostic tests for severe and common bacterial pathogens and antibiotic resistance will (i) save lives by optimizing treatment, (ii) diminish antibiotic resistance by reducing the use of broad spectrum antibiotics and (iii) decrease overall health costs by preventing or shortening hospitalizations.

SUMMARY OF THE INVENTION

[0017] In accordance with the present invention, there is provided sequence from genomic DNA fragments (size of at least 100 base pairs and all described in the sequence listing) selected either by hybridization from genomic libraries or from data banks and which are specific for the detection of commonly encountered bacterial pathogens (i.e. *Escherichia coli*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa*, *Proteus mirabilis*, *Streptococcus pneumoniae*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Enterococcus faecalis*, *Staphylococcus saprophyticus*, *Streptococcus pyogenes*, *Haemophilus influenzae* and *Moraxella catarrhalis*) in clinical specimens. These bacterial species are associated with approximately 90% of urinary tract infections and with a high percentage of other severe infections including septicemia, meningitis, pneumonia, intraabdominal infections, skin infections and many other severe respiratory tract infections. Overall, the above bacterial species may account for up to 80% of bacterial pathogens isolated in routine microbiology laboratories.

[0018] Synthetic oligonucleotides for hybridization (probes) or DNA amplification (primers) were derived from the above species-specific DNA fragments (ranging in sizes from 0.25 to 5.0 kilobase pairs (kbp)) or from selected data bank sequences (GenBank and EMBL). Bacterial species for which some of the oligonucleotide probes and amplification primers were derived from selected data bank sequences are *Escherichia coli*, *Enterococcus faecalis*, *Streptococcus pyogenes* and *Pseudomonas aeruginosa*. The person skilled in the art understands that the important innovation in this invention is the identification of the species-specific DNA fragments selected either from bacterial genomic libraries by hybridization or from data bank sequences. The selection of oligonucleotides from these fragments suitable for diagnostic purposes is also innovative. Specific and ubiquitous oligonucleotides different from the ones tested in the practice are considered as embodiments of the present invention.

[0019] The development of hybridization (with either fragment or oligonucleotide probes) or of DNA amplification protocols for the detection of pathogens from clinical specimens renders possible a very rapid bacterial identification. This will greatly reduce the time currently required for the identification of pathogens in the clinical laboratory

since these technologies can be applied for bacterial detection and identification directly from clinical specimens with minimum pretreatment of any biological specimens to release bacterial DNA. In addition to being 100% specific, probes and amplification primers allow identification of the bacterial species directly from clinical specimens or, alternatively, from an isolated colony. DNA amplification assays have the added advantages of being faster and more sensitive than hybridization assays, since they allow rapid and exponential in vitro replication of the target segment of DNA from the bacterial genome. Universal probes and amplification primers selected from the 16S or 23S rRNA genes highly conserved among bacteria, which permit the detection of any bacterial pathogens, will serve as a procedure to screen out the numerous negative clinical specimens received in diagnostic laboratories. The use of oligonucleotide probes or primers complementary to characterized bacterial genes encoding resistance to antibiotics to identify commonly encountered and clinically important resistance genes is also under the scope of this invention.

DETAILED DESCRIPTION OF THE INVENTION

[0020] Development of Species-Specific DNA Probes

[0021] DNA fragment probes were developed for the following bacterial species: *Escherichia coli*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa*, *Proteus mirabilis*, *Streptococcus pneumoniae*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Staphylococcus saprophyticus*, *Haemophilus influenzae* and *Moraxella catarrhalis*. (For *Enterococcus faecalis* and *Streptococcus pyogenes*, oligonucleotide sequences were exclusively derived from selected data bank sequences). These species-specific fragments were selected from bacterial genomic libraries by hybridization to DNA from a variety of Gram positive and Gram negative bacterial species (Table 5).

[0022] The chromosomal DNA from each bacterial species for which probes were sought was isolated using standard methods. DNA was digested with a frequently cutting restriction enzyme such as Sau3AI and then ligated into the bacterial plasmid vector pGEM3Zf (Promega) linearized by appropriate restriction endonuclease digestion. Recombinant plasmids were then used to transform competent *E. coli* strain DH5 α thereby yielding a genomic library. The plasmid content of the transformed bacterial cells was analyzed using standard methods. DNA fragments of target bacteria ranging in size from 0.25 to 5.0 kilobase pairs (kbp) were cut out from the vector by digestion of the recombinant plasmid with various restriction endonucleases. The insert was separated from the vector by agarose gel electrophoresis and purified in low melting point agarose gels. Each of the purified fragments of bacterial genomic DNA was then used as a probe for specificity tests.

[0023] For each given species, the gel-purified restriction fragments of unknown coding potential were labeled with the radioactive nucleotide α -³²P(dATP) which was incorporated into the DNA fragment by the random priming labeling reaction. Non-radioactive modified nucleotides could also be incorporated into the DNA by this method to serve as a label.

[0024] Each DNA fragment probe (i.e. a segment of bacterial genomic DNA of at least 100 bp in length cut out

from clones randomly selected from the genomic library) was then tested for its specificity by hybridization to DNAs from a variety of bacterial species (Table 5). The double-stranded labeled DNA probe was heat-denatured to yield labeled single-stranded DNA which could then hybridize to any single-stranded target DNA fixed onto a solid support or in solution. The target DNAs consisted of total cellular DNA from an array of bacterial species found in clinical samples (Table 5). Each target DNA was released from the bacterial cells and denatured by conventional methods and then irreversibly fixed onto a solid support (e.g. nylon or nitrocellulose membranes) or free in solution. The fixed single-stranded target DNAs were then hybridized with the single-stranded probe. Pre-hybridization, hybridization and post-hybridization conditions were as follows: (i) Pre-hybridization; in 1 M NaCl+10% dextran sulfate+1% SDS (sodium dodecyl sulfate)+1 μ g/ml salmon sperm DNA at 650° C. for 15 min. (ii) Hybridization; in fresh pre-hybridization solution containing the labeled probe at 650° C. overnight. (iii) Post-hybridization; washes twice in 3 \times SSC containing 1% SDS (1 \times SSC is 0.15M NaCl, 0.015M NaCitrate) and twice in 0.1 \times SSC containing 0.1% SDS; all washes were at 650° C. for 15 min. Autoradiography of washed filters allowed the detection of selectively hybridized probes. Hybridization of the probe to a specific target DNA indicated a high degree of similarity between the nucleotide sequence of these two DNAs. Species-specific DNA fragments selected from various bacterial genomic libraries ranging in size from 0.25 to 5.0 kbp were isolated for 10 common bacterial pathogens (Table 6) based on hybridization to chromosomal DNAs from a variety of bacteria performed as described above. All of the bacterial species tested (66 species listed in Table 5) were likely to be pathogens associated with common infections or potential contaminants which can be isolated from clinical specimens. A DNA fragment probe was considered specific only when it hybridized solely to the pathogen from which it was isolated. DNA fragment probes found to be specific were subsequently tested for their ubiquity (i.e. ubiquitous probes recognized most isolates of the target species) by hybridization to bacterial DNAs from approximately 10 to 80 clinical isolates of the species of interest (Table 6). The DNAs were denatured, fixed onto nylon membranes and hybridized as described above.

[0025] Sequencing of the Species-Specific Fragment Probes

[0026] The nucleotide sequence of the totality or of a portion of the species-specific DNA fragments isolated (Table 6) was determined using the dideoxynucleotide termination sequencing method which was performed using Sequenase (USB Biochemicals) or T7 DNA polymerase (Pharmacia). These nucleotide sequences are shown in the sequence listing. Alternatively, sequences selected from data banks (GenBank and EMBL) were used as sources of oligonucleotides for diagnostic purposes for *Escherichia coli*, *Enterococcus faecalis*, *Streptococcus pyogenes* and *Pseudomonas aeruginosa*. For this strategy, an array of suitable oligonucleotide primers or probes derived from a variety of genomic DNA fragments (size of more than 100 bp) selected from data banks was tested for their specificity and ubiquity in PCR and hybridization assays as described later. It is important to note that the data bank sequences were selected based on their potential of being species-specific according to available sequence information. Only

data bank sequences from which species-specific oligonucleotides could be derived are included in this invention.

[0027] Oligonucleotide probes and amplification primers derived from species-specific fragments selected from the genomic libraries or from data bank sequences were synthesized using an automated DNA synthesizer (Millipore). Prior to synthesis, all oligonucleotides (probes for hybridization and primers for DNA amplification) were evaluated for their suitability for hybridization or DNA amplification by polymerase chain reaction (PCR) by computer analysis using standard programs (e.g. Genetics Computer Group (GCG) and Oligo™ 4.0 (National Biosciences)). The potential suitability of the PCR primer pairs was also evaluated prior to the synthesis by verifying the absence of unwanted features such as long stretches of one nucleotide, a high proportion of G or C residues at the 3' end and a 3'-terminal T residue (Persing et al, 1993. *Diagnostic Molecular Microbiology: Principles and Applications*, American Society for Microbiology, Washington, D.C.).

[0028] Hybridization with Oligonucleotide Probes

[0029] In hybridization experiments, oligonucleotides (size less than 100 nucleotides) have some advantages over DNA fragment probes for the detection of bacteria such as ease of preparation in large quantities, consistency in results from batch to batch and chemical stability. Briefly, for the hybridizations, oligonucleotides were 5' end-labeled with the radionucleotide $\gamma^{32}\text{P}(\text{ATP})$ using T4 polynucleotide kinase (Pharmacia). The unincorporated radionucleotide was removed by passing the labeled single-stranded oligonucleotide through a Sephadex G50 column. Alternatively, oligonucleotides were labeled with biotin, either enzymatically at their 3' ends or incorporated directly during synthesis at their 5' ends, or with digoxigenin. It will be appreciated by the person skilled in the art that labeling means other than the three above labels may be used.

[0030] The target DNA was denatured, fixed onto a solid support and hybridized as previously described for the DNA fragment probes. Conditions for pre-hybridization and hybridization were as described earlier. Post-hybridization washing conditions were as follows: twice in 3×SSC containing 1% SDS, twice in 2×SSC containing 1% SDS and twice in 1×SSC containing 1% SDS (all of these washes were at 65° C. for 15 min), and a final wash in 0.1×SSC containing 1% SDS at 25° C. for 15 min. For probes labeled with radioactive labels the detection of hybrids was by autoradiography as described earlier. For non-radioactive labels detection may be calorimetric or by chemiluminescence.

[0031] The oligonucleotide probes may be derived from either strand of the duplex DNA. The probes may consist of the bases A, G, C, or T or analogs. The probes may be of any suitable length and may be selected anywhere within the species-specific genomic DNA fragments selected from the genomic libraries or from data bank sequences.

[0032] DNA Amplification

[0033] For DNA amplification by the widely used PCR (polymerase chain reaction) method, primer pairs were derived either from the sequenced species-specific DNA fragments or from data bank sequences or, alternatively, were shortened versions of oligonucleotide probes. Prior to synthesis, the potential primer pairs were analyzed by using

the program oligo™ 4.0 (National Biosciences) to verify that they are likely candidates for PCR amplifications.

[0034] During DNA amplification by PCR, two oligonucleotide primers binding respectively to each strand of the denatured double-stranded target DNA from the bacterial genome are used to amplify exponentially in vitro the target DNA by successive thermal cycles allowing denaturation of the DNA, annealing of the primers and synthesis of new targets at each cycle (Persing et al, 1993. *Diagnostic Molecular Microbiology: Principles and Applications*, American Society for Microbiology, Washington, D.C.). Briefly, the PCR protocols were as follows. Clinical specimens or bacterial colonies were added directly to the 50 μL PCR reaction mixtures containing 50 mM KCl, 10 mM Tris-HCl pH 8.3, 2.5 mM MgCl_2 , 0.4 μM of each of the two primers, 200 μM of each of the four dNTPs and 1.25 Units of Taq DNA polymerase (Perkin Elmer). PCR reactions were then subjected to thermal cycling (3 min at 95° C. followed by 30 cycles of 1 second at 95° C. and 1 second at 55° C.) using a Perkin Elmer 480™ thermal cycler and subsequently analyzed by standard ethidium bromide-stained agarose gel electrophoresis. It is clear that other methods for the detection of specific amplification products, which may be faster and more practical for routine diagnosis, may be used. Such methods may be based on the detection of fluorescence after amplification (e.g. TaqMan™ system from Perkin Elmer or Amplisensor™ from Biotronics) or liquid hybridization with an oligonucleotide probe binding to internal sequences of the specific amplification product. These novel probes can be generated from our species-specific fragment probes. Methods based on the detection of fluorescence are particularly promising for utilization in routine diagnosis as they are, very rapid and quantitative and can be automated.

[0035] To assure PCR efficiency, glycerol or dimethyl sulfoxide (DMSO) or other related solvents, can be used to increase the sensitivity of the PCR and to overcome problems associated with the amplification of target with a high GC content or with strong secondary structures. The concentration ranges for glycerol and DMSO are 5-15% (v/v) and 3-10% (v/v), respectively. For the PCR reaction mixture, the concentration ranges for the amplification primers and the MgCl_2 are 0.1-1.0 μM and 1.5-3.5 mM, respectively. Modifications of the standard PCR protocol using external and nested primers (i.e. nested PCR) or using more than one primer pair (i.e. multiplex PCR) may also be used (Persing et al, 1993. *Diagnostic Molecular Microbiology: Principles and Applications*, American Society for Microbiology, Washington, D.C.). For more details about the PCR protocols and amplicon detection methods see examples 7 and 8.

[0036] The person skilled in the art of DNA amplification knows the existence of other rapid amplification procedures such as ligase chain reaction (LCR), transcription-based amplification systems (TAS), self-sustained sequence replication (3SR), nucleic acid sequence-based amplification (NASBA), strand displacement amplification (SDA) and branched DNA (bDNA) (Persing et al, 1993. *Diagnostic Molecular Microbiology: Principles and Applications*, American Society for Microbiology, Washington, D.C.). The scope of this invention is not limited to the use of amplification by PCR, but rather includes the use of any rapid nucleic acid amplification methods or any other procedures which may be used to increase rapidity and sensitivity of the

tests. Any oligonucleotides suitable for the amplification of nucleic acid by approaches other than PCR and derived from the species-specific fragments and from selected antibiotic resistance gene sequences included in this document are also under the scope of this invention.

[0037] Specificity and Ubiquity Tests for Oligonucleotide Probes and Primers

[0038] The specificity of oligonucleotide probes, derived either from the sequenced species-specific fragments or from data bank sequences, was tested by hybridization to DNAs from the array of bacterial species listed in Table 5 as previously described. Oligonucleotides found to be specific were subsequently tested for their ubiquity by hybridization to bacterial DNAs from approximately 80 isolates of the target species as described for fragment probes. Probes were considered ubiquitous when they hybridized specifically with the DNA from at least 80% of the isolates. Results for specificity and ubiquity tests with the oligonucleotide probes are summarized in Table 6. The specificity and ubiquity of the amplification primer pairs were tested directly from cultures (see example 7) of the same bacterial strains. For specificity and ubiquity tests, PCR assays were performed directly from bacterial colonies of approximately 80 isolates of the target species. Results are summarized in Table 7. All specific and ubiquitous oligonucleotide probes and amplification primers for each of the 12 bacterial species investigated are listed in Annexes I and II, respectively. Divergence in the sequenced DNA fragments can occur and, insofar as the divergence of these sequences or a part thereof does not affect the specificity of the probes or amplification primers, variant bacterial DNA is under the scope of this invention.

[0039] Universal Bacterial Detection

[0040] In the routine microbiology laboratory a high percentage of clinical specimens sent for bacterial identification is negative (Table 4). For example, over a 2 year period, around 80% of urine specimens received by the laboratory at the "Centre Hospitalier de l' Université Laval (CHUL)" were negative (i.e. $<10^7$ CFU/L) (Table 3). Testing clinical samples with universal probes or universal amplification primers to detect the presence of bacteria prior to specific identification and screen out the numerous negative specimens is thus useful as it saves costs and may rapidly orient the clinical management of the patients. Several oligonucleotides and amplification primers were therefore synthesized from highly conserved portions of bacterial 16S or 23S ribosomal RNA gene sequences available in data banks (Annexes III and IV). In hybridization tests, a pool of seven oligonucleotides (Annex I; Table 6) hybridized strongly to DNA from all bacterial species listed in Table 5. This pool of universal probes labeled with radionucleotides or with any other modified nucleotides is consequently very useful for detection of bacteria in urine samples with a sensitivity range of $\geq 10^7$ CFU/L. These probes can also be applied for bacterial detection in other clinical samples.

[0041] Amplification primers also derived from the sequence of highly conserved ribosomal RNA genes were used as an alternative strategy for universal bacterial detection directly from clinical specimens (Annex IV; Table 7). The DNA amplification strategy was developed to increase the sensitivity and the rapidity of the test. This amplification test was ubiquitous since it specifically amplified DNA from 23 different bacterial species encountered in clinical specimens.

[0042] Well-conserved bacterial genes other than ribosomal RNA genes could also be good candidates for universal bacterial detection directly from clinical specimens. Such genes may be associated with processes essential for bacterial survival (e.g. protein synthesis, DNA synthesis, cell division or DNA repair) and could therefore be highly conserved during evolution. We are working on these candidate genes to develop new rapid tests for the universal detection of bacteria directly from clinical specimens.

[0043] Antibiotic Resistance Genes

[0044] Antimicrobial resistance complicates treatment and often leads to therapeutic failures. Furthermore, overuse of antibiotics inevitably leads to the emergence of bacterial resistance. Our goal is to provide the clinicians, within one hour, the needed information to prescribe optimal treatments. Besides the rapid identification of negative clinical specimens with DNA-based tests for universal bacterial detection and the identification of the presence of a specific pathogen in the positive specimens with DNA-based tests for specific bacterial detection, the clinicians also need timely information about the ability of the bacterial pathogen to resist antibiotic treatments. We feel that the most efficient strategy to evaluate rapidly bacterial resistance to antimicrobials is to detect directly from the clinical specimens the most common and important antibiotic resistance genes (i.e. DNA-based tests for the detection of antibiotic resistance genes). Since the sequence from the most important and common bacterial antibiotic resistance genes are available from data banks, our strategy is to use the sequence from a portion or from the entire gene to design specific oligonucleotides which will be used as a basis for the development of rapid DNA-based tests. The sequence from the bacterial antibiotic resistance genes selected on the basis of their clinical relevance (i.e. high incidence and importance) is given in the sequence listing. Table 8 summarizes some characteristics of the selected antibiotic resistance genes.

EXAMPLES

[0045] The following examples are intended to be illustrative of the various methods and compounds of the invention.

Example 1

[0046] Isolation and cloning of fragments. Genomic DNAs from *Escherichia coli* strain ATCC 25922, *Klebsiella pneumoniae* strain CK2, *Pseudomonas aeruginosa* strain ATCC 27853, *Proteus mirabilis* strain ATCC 35657, *Streptococcus pneumoniae* strain ATCC 27336, *Staphylococcus aureus* strain ATCC 25923, *Staphylococcus epidermidis* strain ATCC 12228, *Staphylococcus saprophyticus* strain ATCC 15305, *Haemophilus influenzae* reference strain Rd and *Moraxella catarrhalis* strain ATCC 53879 were prepared using standard procedures. It is understood that the bacterial genomic DNA may have been isolated from strains other than the ones mentioned above. (For *Enterococcus faecalis* and *Streptococcus pyogenes* oligonucleotide sequences were derived exclusively from data banks). Each DNA was digested with a restriction enzyme which frequently cuts DNA such as Sau3AI. The resulting DNA fragments were ligated into a plasmid vector (pGEM3Zf) to create recombinant plasmids and transformed into compe-

tent *E. coli* cells (DH5 α). It is understood that the vectors and corresponding competent cells should not be limited to the ones herein above specifically exemplified. The objective of obtaining recombinant plasmids and transformed cells is to provide an easily reproducible source of DNA fragments useful as probes. Therefore, insofar as the inserted fragments are specific and selective for the target bacterial DNA, any recombinant plasmids and corresponding transformed host cells are under the scope of this invention. The plasmid content of the transformed bacterial cells was analyzed using standard methods. DNA fragments from target bacteria ranging in size from 0.25 to 5.0 kbp were cut out from the vector by digestion of the recombinant plasmid with various restriction endonucleases. The insert was separated from the vector by agarose gel electrophoresis and purified in a low melting point agarose gel. Each of the purified fragments was then used for specificity tests.

[0047] Labeling of DNA fragment probes. The label used was α -³²P(dATP), a radioactive nucleotide which can be incorporated enzymatically into a double-stranded DNA molecule. The fragment of interest is first denatured by heating at 95° C. for 5 min, then a mixture of random primers is allowed to anneal to the strands of the fragments. These primers, once annealed, provide a starting point for synthesis of DNA. DNA polymerase, usually the Klenow fragment, is provided along with the four nucleotides, one of which is radioactive. When the reaction is terminated, the mixture of new DNA molecules is once again denatured to provide radioactive single-stranded DNA molecules (i.e. the probe). As mentioned earlier, other modified nucleotides may be used to label the probes.

[0048] Specificity and ubiquity tests for the DNA fragment probes. Species-specific DNA fragments ranging in size from 0.25 to 5.0 kbp were isolated for 10 common bacterial pathogens (Table 6) based on hybridization to chromosomal DNAs from a variety of bacteria. Samples of whole cell DNA for each bacterial strain listed in Table 5 were transferred onto a nylon membrane using a dot blot apparatus, washed and denatured before being irreversibly fixed. Hybridization conditions were as described earlier. A DNA fragment probe was considered specific only when it hybridized solely to the pathogen from which it was isolated. Labeled DNA fragments hybridizing specifically only to target bacterial species (i.e. specific) were then tested for their ubiquity by hybridization to DNAs from approximately 10 to 80 isolates of the species of interest as described earlier. The conditions for pre-hybridization, hybridization and post-hybridization washes were as described earlier. After autoradiography (or other detection means appropriate for the non-radioactive label used), the specificity of each individual probe can be determined. Each probe found to be specific (i.e. hybridizing only to the DNA from the bacterial species from which it was isolated) and ubiquitous (i.e. hybridizing to most isolates of the target species) was kept for further experimentations.

Example 2

[0049] Same as example 1 except that testing of the strains is by colony hybridization. The bacterial strains were inoculated onto a nylon membrane placed on nutrient agar. The membranes were incubated at 37° C. for two hours and then bacterial lysis and DNA denaturation were carried out according to standard procedures. DNA hybridization was performed as described earlier.

Example 3

[0050] Same as example 1 except that bacteria were detected directly from clinical samples. Any biological samples were loaded directly onto a dot blot apparatus and cells were lysed in situ for bacterial detection. Blood samples should be heparinized in order to avoid coagulation interfering with their convenient loading on a dot blot apparatus.

Example 4

[0051] Nucleotide sequencing of DNA fragments. The nucleotide sequence of the totality or a portion of each fragment found to be specific and ubiquitous (Example 1) was determined using the dideoxynucleotide termination sequencing method (Sanger et al., 1977, Proc. Natl. Acad. Sci. USA. 74:5463-5467). These DNA sequences are shown in the sequence listing. Oligonucleotide probes and amplification primers were selected from these nucleotide sequences, or alternatively, from selected data banks and were then synthesized on an automated Biosearch synthesizer (Millipore™) using phosphoramidite chemistry.

[0052] Labeling of oligonucleotides. Each oligonucleotide was 5' end-labeled with γ -³²P-ATP by the T4 polynucleotide kinase (Pharmacia) as described earlier. The label could also be non-radioactive.

[0053] Specificity test for oligonucleotide probes. All labeled oligonucleotide probes were tested for their specificity by hybridization to DNAs from a variety of Gram positive and Gram negative bacterial species as described earlier. Species-specific probes were those hybridizing only to DNA from the bacterial species from which it was isolated. Oligonucleotide probes found to be specific were submitted to ubiquity tests as follows.

[0054] Ubiquity test for oligonucleotide probes. Specific oligonucleotide probes were then used in ubiquity tests with approximately 80 strains of the target species. Chromosomal DNAs from the isolates were transferred onto nylon membranes and hybridized with labeled oligonucleotide probes as described for specificity tests. The batteries of approximately 80 isolates constructed for each target species contain reference ATCC strains as well as a variety of clinical isolates obtained from various sources. Ubiquitous probes were those hybridizing to at least 80% of DNAs from the battery of clinical isolates of the target species. Examples of specific and ubiquitous oligonucleotide probes are listed in Annex 1.

Example 5

[0055] Same as example 4 except that a pool of specific oligonucleotide probes is used for bacterial identification (i) to increase sensitivity and assure 100% ubiquity or (ii) to identify simultaneously more than one bacterial species. Bacterial identification could be done from isolated colonies or directly from clinical specimens.

Example 6

[0056] PCR amplification. The technique of PCR was used to increase sensitivity and rapidity of the tests. The PCR primers used were often shorter derivatives of the extensive sets of oligonucleotides previously developed for

hybridization assays (Table 6). The sets of primers were tested in PCR assays performed directly from a bacterial colony or from a bacterial suspension (see Example 7) to determine their specificity and ubiquity (Table 7). Examples of specific and ubiquitous PCR primer pairs are listed in annex II.

[0057] Specificity and ubiquity tests for amplification primers. The specificity of all selected PCR primer pairs was tested against the battery of Gram negative and Gram positive bacteria used to test the oligonucleotide probes (Table 5). Primer pairs found specific for each species were then tested for their ubiquity to ensure that each set of primers could amplify at least 80% of DNAs from a battery of approximately 80 isolates of the target species. The batteries of isolates constructed for each species contain reference ATCC strains and various clinical isolates representative of the clinical diversity for each species.

[0058] Standard precautions to avoid false positive PCR results should be taken. Methods to inactivate PCR amplification products such as the inactivation by uracil-N-glycosylase may be used to control PCR carryover.

Example 7

[0059] Amplification directly from a bacterial colony or suspension. PCR assays were performed either directly from a bacterial colony or from a bacterial suspension, the latter being adjusted to a standard McFarland 0.5 (corresponds to 1.5×10^8 bacteria/mL). In the case of direct amplification from a colony, a portion of the colony was transferred directly to a 50 μ L PCR reaction mixture (containing 50 mM KCl, 10 mM Tris pH 8.3, 2.5 mM MgCl₂, 0.4 μ M of each of the two primers, 200 μ M of each of the four dNTPs and 1.25 Unit of Taq DNA polymerase (Perkin Elmer)) using a plastic rod. For the bacterial suspension, 4 μ L of the cell suspension was added to 46 μ L of the same PCR reaction mixture. For both strategies, the reaction mixture was overlaid with 50 μ L of mineral oil and PCR amplifications were carried out using an initial denaturation step of 3 min. at 95° C. followed by 30 cycles consisting of a 1 second denaturation step at 95° C. and of a 1 second annealing step at 55° C. in a Perkin Elmer 480™ thermal cycler. PCR amplification products were then analyzed by standard agarose gel (2%) electrophoresis. Amplification products were visualized in agarose gels containing 2.5 μ g/mL of ethidium bromide under UV at 254 nm. The entire PCR assay can be completed in approximately one hour.

[0060] Alternatively, amplification from bacterial cultures was performed as described above but using a “hot start” protocol. In that case, an initial reaction mixture containing the target DNA, primers and dNTPs was heated at 85° C. prior to the addition of the other components of the PCR reaction mixture. The final concentration of all reagents was as described above. Subsequently, the PCR reactions were submitted to thermal cycling and analysis as described above.

Example 8

[0061] Amplification directly from clinical specimens. For amplification from urine specimens, 4 μ L of undiluted or diluted (1:10) urine was added directly to 46 μ L of the above PCR reaction mixture and amplified as described earlier.

[0062] To improve bacterial cell lysis and eliminate the PCR inhibitory effects of clinical specimens, samples were routinely diluted in lysis buffer containing detergent(s). Subsequently, the lysate was added directly to the PCR reaction mixture. Heat treatments of the lysates, prior to DNA amplification, using the thermocycler or a microwave oven could also be performed to increase the efficiency of cell lysis.

[0063] Our strategy is to develop rapid and simple protocols to eliminate PCR inhibitory effects of clinical specimens and lyse bacterial cells to perform DNA amplification directly from a variety of biological samples. PCR has the advantage of being compatible with crude DNA preparations. For example, blood, cerebrospinal fluid and sera may be used directly in PCR assays after a brief heat treatment. We intend to use such rapid and simple strategies to develop fast protocols for DNA amplification from a variety of clinical specimens.

Example 9

[0064] Detection of antibiotic resistance genes. The presence of specific antibiotic resistance genes which are frequently encountered and clinically relevant is identified using the PCR amplification or hybridization protocols described in previous sections. Specific oligonucleotides used as a basis for the DNA-based tests are selected from the antibiotic resistance gene sequences. These tests can be performed either directly from clinical specimens or from a bacterial colony and should complement diagnostic tests for specific bacterial identification.

Example 10

[0065] Same as examples 7 and 8 except that assays were performed by multiplex PCR (i.e. using several pairs of primers in a single PCR reaction) to (i) reach an ubiquity of 100% for the specific target pathogen or (ii) to detect simultaneously several species of bacterial pathogens.

[0066] For example, the detection of *Escherichia coli* requires three pairs of PCR primers to assure a ubiquity of 100%. Therefore, a multiplex PCR assay (using the “hot-start” protocol (Example 7)) with those three primer pairs was developed. This strategy was also used for the other bacterial pathogens for which more than one primer pair was required to reach an ubiquity of 100%.

[0067] Multiplex PCR assays could also be used to (i) detect simultaneously several bacterial species or, alternatively, (ii) to simultaneously identify the bacterial pathogen and detect specific antibiotic resistance genes either directly from a clinical specimen or from a bacterial colony.

[0068] For these applications, amplicon detection methods should be adapted to differentiate the various amplicons produced. Standard agarose gel electrophoresis could be used because it discriminates the amplicons based on their sizes. Another useful strategy for this purpose would be detection using a variety of fluorochromes emitting at different wavelengths which are each coupled with a specific oligonucleotide linked to a fluorescence quencher which is degraded during amplification to release the fluorochrome (e.g. TaqMan™, Perkin Elmer).

Example 11

[0069] Detection of amplification Products. The person skilled in the art will appreciate that alternatives other than

standard agarose gel electrophoresis (Example 7) may be used for the revelation of amplification products. Such methods may be based on the detection of fluorescence after amplification (e.g. Amplisensor™, Biotronics; TaqMan™) or other labels such as biotin (SHARP Signal™ system, Digene Diagnostics). These methods are quantitative and easily automated. One of the amplification primers or an internal oligonucleotide probe specific to the amplicon(s) derived from the species-specific fragment probes is coupled with the fluorochrome or with any other label. Methods based on the detection of fluorescence are particularly suitable for diagnostic tests since they are rapid and flexible as fluorochromes emitting different wavelengths are available (Perkin Elmer).

Example 12

[0070] Species-specific, universal and antibiotic resistance gene amplification primers can be used in other rapid amplification procedures such as the ligase chain reaction (LCR), transcription-based amplification systems (TAS), self-sustained sequence replication (3SR), nucleic acid sequence-based amplification (NASBA), strand displacement amplification (SDA) and branched DNA (bDNA) or any other methods to increase the sensitivity of the test. Amplifications can be performed from an isolated bacterial colony or directly from clinical specimens. The scope of this invention is therefore not limited to the use of PCR but rather includes the use of any procedures to specifically identify bacterial DNA and which may be used to increase rapidity and sensitivity of the tests.

Example 13

[0071] A test kit would contain sets of probes specific for each bacterium as well as a set of universal probes. The kit is provided in the form of test components, consisting of the set of universal probes labeled with non-radioactive labels as well as labeled specific probes for the detection of each bacterium of interest in specific clinical samples. The kit will also include test reagents necessary to perform the pre-hybridization, hybridization, washing steps and hybrid detection. Finally, test components for the detection of known antibiotic resistance genes (or derivatives therefrom) will be included. Of course, the kit will include standard samples to be used as negative and positive controls for each hybridization test.

[0072] Components to be included in the kits will be adapted to each specimen type and to detect pathogens commonly encountered in that type of specimen. Reagents for the universal detection of bacteria will also be included. Based on the sites of infection, the following kits for the specific detection of pathogens may be developed:

[0073] A kit for the universal detection of bacterial pathogens from most clinical specimens which contains sets of probes specific for highly conserved regions of the bacterial genomes.

[0074] A kit for the detection of bacterial pathogens retrieved from urine samples, which contains eight specific test components (sets of probes for the

detection of *Escherichia coli*, *Enterococcus faecalis*, *Klebsiella pneumoniae*, *Proteus mirabilis*, *Pseudomonas aeruginosa*, *Staphylococcus saprophyticus*, *Staphylococcus aureus* and *Staphylococcus epidermidis*).

[0075] A kit for the detection of respiratory pathogens which contains seven specific test components (sets of probes for detecting *Streptococcus pneumoniae*, *Moraxella catarrhalis*, *Haemophilus influenzae*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa*, *Streptococcus pyogenes* and *Staphylococcus aureus*).

[0076] A kit for the detection of pathogens retrieved from blood samples, which contains eleven specific test components (sets of probes for the detection of *Streptococcus pneumoniae*, *Moraxella catarrhalis*, *Haemophilus influenzae*, *Proteus mirabilis*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa*, *Escherichia coli*, *Enterococcus faecalis*, *Staphylococcus aureus*, *Streptococcus pyogenes* and *Staphylococcus epidermidis*).

[0077] A kit for the detection of pathogens causing meningitis, which contains four specific test components (sets of probes for the detection of *Haemophilus influenzae*, *Streptococcus pneumoniae*, *Escherichia coli* and *Pseudomonas aeruginosa*).

[0078] A kit for the detection of clinically important antibiotic resistance genes which contains sets of probes for the specific detection of at least one of the 19 following genes associated with bacterial resistance : bla_{tem}, bla_{rob}, bla_{shv}, aadB, aacC1, aacC2, aacC3, aacA4, mecA, vanA, vanH, vanX, satA, aacA-aphD, vat, vga, msrA, sul and int.

[0079] Other kits adapted for the detection of pathogens from skin, abdominal wound or any other clinically relevant kits will be developed.

Example 14

[0080] Same as example 13 except that the test kits contain all reagents and controls to perform DNA amplification assays. Diagnostic kits will be adapted for amplification by PCR (or other amplification methods) performed directly either from clinical specimens or from a bacterial colony. Components required for universal bacterial detection, bacterial identification and antibiotic resistance genes detection will be included.

[0081] Amplification assays could be performed either in tubes or in microtitration plates having multiple wells. For assays in plates, the wells will be coated with the specific amplification primers and control DNAs and the detection of amplification products will be automated. Reagents and amplification primers for universal bacterial detection will be included in kits for tests performed directly from clinical specimens. Components required for bacterial identification and antibiotic resistance gene detection will be included in

kits for testing directly from colonies as well as in kits for testing directly from clinical specimens.

[0082] The kits will be adapted for use with each type of specimen as described in example 13 for hybridization-based diagnostic kits.

Example 15

[0083] It is understood that the use of the probes and amplification primers described in this invention for bacterial detection and identification is not limited to clinical microbiology applications. In fact, we feel that other sectors could also benefit from these new technologies. For example, these tests could be used by industries for quality control of food, water, pharmaceutical products or other products requiring microbiological control. These tests could also be applied to detect and identify bacteria in biological samples from organisms other than humans (e.g. other primates, mammals, farm animals and live stocks). These diagnostic tools could also be very useful for research purposes including clinical trials and epidemiological studies.

TABLE 1

Distribution of urinary isolates from positive urine samples ($\geq 10^7$ CFU/L) at the Centre Hospitalier de l'Université Laval (CHUL) for the 1992–1994 period.				
Organisms	% of isolates			
	Nov 92 n = 267 ^a	April 93 n = 265	July 93 n = 238	Jan 94 n = 281
<i>Escherichia coli</i>	53.2	51.7	53.8	54.1
<i>Enterococcus faecalis</i>	13.8	12.4	11.7	11.4
<i>Klebsiella pneumoniae</i>	6.4	6.4	5.5	5.3
<i>Staphylococcus epidermidis</i>	7.1	7.9	3.0	6.4
<i>Proteus mirabilis</i>	2.6	3.4	3.8	2.5
<i>Pseudomonas aeruginosa</i>	3.7	3.0	5.0	2.9
<i>Staphylococcus saprophyticus</i>	3.0	1.9	5.4	1.4
Others ^b	10.2	13.3	11.8	16.0

^an = total number of isolates for the indicated month.

^bSee table 2.

[0084]

TABLE 2

Distribution of uncommon ^a urinary isolates from positive urine samples ($\geq 10^7$ CFU/L) at the Centre Hospitalier de l'Université Laval (CHUL) for the 1992–1994 period.				
Organisms ^a	% of isolates			
	Nov 92	April 93	July 93	Jan 94
<i>Staphylococcus aureus</i>	0.4	1.1	1.3	1.4
<i>Staphylococcus</i> spp.	2.2	4.9	1.7	6.0
<i>Micrococcus</i> spp.	0.0	0.0	0.4	0.7
<i>Enterococcus faecium</i>	0.4	0.4	1.3	1.4
<i>Citrobacter</i> spp.	1.4	0.8	0.4	0.7
<i>Enterobacter</i> spp.	1.5	1.1	1.3	1.4
<i>Klebsiella oxytoca</i>	1.1	1.5	2.5	1.8
<i>Serratia</i> spp.	0.8	0.0	0.5	0.0
<i>Proteus</i> spp.	0.4	0.4	0.0	1.1

TABLE 2-continued

Distribution of uncommon ^a urinary isolates from positive urine samples ($\geq 10^7$ CFU/L) at the Centre Hospitalier de l'Université Laval (CHUL) for the 1992–1994 period.				
Organisms ^a	% of isolates			
	Nov 92	April 93	July 93	Jan 94
<i>Morganella</i> and <i>Providencia</i>	0.4	0.8	0.4	0.0
<i>Hafnia alvei</i>	0.8	0.0	0.0	0.0
NFB ^b (<i>Stenotrophomonas</i> , <i>Acinetobacter</i>)	0.0	0.4	1.3	1.1
<i>Candida</i> spp.	0.8	1.9	0.7	0.4

^aUncommon urinary isolates are those identified as "Others" in Table 1.

^bNFB: non fermentative bacilli.

[0085]

TABLE 3

Distribution of positive ^a (bacterial count $\geq 10^7$ CFU/L) and negative (bacterial count $< 10^7$ CFU/L) urine specimens tested at the Centre Hospitalier de l'Université Laval (CHUL) for the 1992–1994 period.				
Specimens	Number of isolates (%)			
	Nov 92	April 93	July 93	Jan 94
received:	1383(100)	1338(100)	1139(100)	1345(100)
positive:	267(19.3)	265(19.8)	238(20.9)	281(20.9)
negative:	1116(80.7)	1073(80.2)	901(79.1)	1064(79.1)

^aBased on standard diagnostic methods, the minimal number of bacterial pathogens in urine samples to indicate an urinary tract infection is normally 10^7 CFU/L.

[0086]

TABLE 4

Distribution of positive and negative clinical specimens tested in the Microbiology Laboratory of the CHUL.			
Clinical specimens ^a	No. of samples tested	% of positive specimens	% of negative specimens
Urine	17,981	19.4	80.6
Haemoculture/marrow	10,010	6.9	93.1
Sputum	1,266	68.4	31.6
Superficial pus	1,136	72.3	27.7
Cerebrospinal fluid	553	1.0	99.0
Synovial fluid-articular	523	2.7	97.3
Bronch./Trach./Amyg./Throat	502	56.6	43.4
Deep pus	473	56.8	43.2
Ears	289	47.1	52.9
Pleural and pericardial fluid	132	1.0	99.0
Peritoneal fluid	101	28.6	71.4

^aSpecimens tested from February 1994 to January 1995.

[0087]

TABLE 5

Bacterial species (66) used for testing the specificity of DNA fragment probes, oligonucleotide probes and PCR primers.			
Bacterial species	Number of strains tested	Bacterial species	Number of strains tested
Gram negative:		Gram positive:	
<i>Proteus mirabilis</i>	5	<i>Streptococcus pneumoniae</i>	7
<i>Klebsiella pneumoniae</i>	5	<i>Streptococcus salivarius</i>	2
<i>Pseudomonas aeruginosa</i>	5	<i>Streptococcus viridans</i>	2
<i>Escherichia coli</i>	5	<i>Streptococcus pyogenes</i>	2
<i>Moraxella catarrhalis</i>	5	<i>Staphylococcus aureus</i>	2
<i>Proteus vulgaris</i>	2	<i>Staphylococcus epidermidis</i>	2
<i>Morganella morganii</i>	2	<i>Staphylococcus saprophyticus</i>	5
<i>Enterobacter cloacae</i>	2	Micrococcus species	2
<i>Providencia stuartii</i>	1	Corynebacterium species	2
Providencia species	1	<i>Streptococcus</i> groupe B	2
<i>Enterobacter agglomerans</i>	2	<i>Staphylococcus simulans</i>	2
<i>Providencia rettgeri</i>	2	<i>Staphylococcus ludgunensis</i>	1
<i>Neisseria mucosa</i>	1	<i>Staphylococcus capitis</i>	2
<i>Providencia alcalifaciens</i>	1	<i>Staphylococcus haemolyticus</i>	2
<i>Providencia rustigianii</i>	1	<i>Staphylococcus hominis</i>	2
<i>Burkholderia cepacia</i>	2	<i>Enterococcus faecalis</i>	2
<i>Enterobacter aerogenes</i>	2	<i>Enterococcus faecium</i>	1
<i>Stenotrophomonas maltophilia</i>	2	<i>Staphylococcus warneri</i>	1
<i>Pseudomonas fluorescens</i>	1	<i>Enterococcus durans</i>	1
<i>Comamonas acidovorans</i>	2	<i>Streptococcus bovis</i>	1
<i>Pseudomonas putida</i>	2	Diptheroids	2
<i>Haemophilus influenzae</i>	5	<i>Lactobacillus acidophilus</i>	1
<i>Haemophilus parainfluenzae</i>	2		
<i>Bordetella pertussis</i>	2		
<i>Haemophilus parahaemolyticus</i>	2		
<i>Haemophilus haemolyticus</i>	2		
<i>Haemophilus aegyptius</i>	1		
<i>Kingella indologenes</i>	1		
<i>Moraxella atlantae</i>	1		
<i>Neisseria caviae</i>	1		
<i>Neisseria subflava</i>	1		
<i>Moraxella urethralis</i>	1		
<i>Shigella sonnei</i>	1		
<i>Shigella flexneri</i>	1		
<i>Klebsiella oxytoca</i>	2		
<i>Serratia marcescens</i>	2		
<i>Salmonella typhimurium</i>	1		
<i>Yersinia enterocolitica</i>	1		
<i>Acinetobacter calcoaceticus</i>	1		
<i>Acinetobacter lwoffii</i>	1		
<i>Haftnia alvei</i>	2		
<i>Citrobacter diversus</i>	1		
<i>Citrobacter freundii</i>	1		
Salmonella species	1		

[0088]

TABLE 6

Organisms ^a	Number of fragment probes ^b			Number of oligonucleotide probes		
	Tested	Specific	Ubiquitous ^c	Synthesized	Specific	Ubiquitous ^c
<i>E. coli</i> ^d	—	—	—	20	12	9 ^f
<i>E. coli</i>	14	2	2 ^e	—	—	—
<i>K. pneumoniae</i> ^d	—	—	—	15	1	1
<i>K. pneumoniae</i>	33	3	3	18	12	8
<i>P. mirabilis</i> ^d	—	—	—	3	3	2
<i>P. mirabilis</i>	14	3	3 ^e	15	8	7
<i>P. aeruginosa</i> ^d	—	—	—	26	13	9
<i>P. aeruginosa</i>	6	2	2 ^e	6	0	0

TABLE 6-continued

Organisms ^a	Number of fragment probes ^b			Number of oligonucleotide probes		
	Tested	Specific	Ubiquitous ^c	Synthesized	Specific	Ubiquitous ^c
<i>S. saprophyticus</i>	7	4	4	20	9	7
<i>H. influenzae</i> ^d	—	—	—	16	2	2
<i>H. influenzae</i>	1	1	1	20	1	1
<i>S. pneumoniae</i> ^d	—	—	—	6	1	1
<i>S. pneumoniae</i>	19	2	2	4	1	1
<i>M. catarrhalis</i>	2	2	2	9	8	8
<i>S. epidermidis</i>	62	1	1	—	—	—
<i>S. aureus</i>	30	1	1	—	—	—
Universal probes ^d	—	—	—	7	—	7 ^B

^aNo DNA fragment or oligonucleotide probes were tested for *E. faecalis* and *S. pyogenes*.

^bSizes of DNA fragments range from 0.25 to 5.0 kbp.

^cA specific probe was considered ubiquitous when at least 80% of isolates of the target species (approximately 80 isolates) were recognized by each specific probe. When 2 or more probes are combined, 100% of the isolates are recognized.

^dThese sequences were selected from data banks.

^eUbiquity tested with approximately 10 isolates of the target species.

^fA majority of probes (8/9) do not discriminate *E. coli* and *Shigella* spp.

^gUbiquity tests with a pool of the 7 probes detected all 66 bacterial species listed in Table 5.

[0089]

TABLE 7

Organism	Primer pair ^a # (SEQ ID NO)	Amplicon size (bp)	Ubiquity ^b	DNA amplification	
				from colonies ^c	from specimens ^d
<i>E. coli</i>	1 ^e (55–56)	107	75/80	+	+
	2 ^e (46–47)	297	77/80	+	+
	3 (42–43)	102	78/80	+	+
	4 (131–132)	134	73/80	+	+
	1 + 3 + 4	—	80/80	+	+
<i>E. faecalis</i>	1 ^e (38–39)	200	71/80	+	+
	2 ^e (40–41)	121	79/80	+	+
	1 + 2	—	80/80	+	+
	1 (67–68)	198	76/80	+	+
<i>K. pneumoniae</i>	2 (61–62)	143	67/80	+	+
	3 ^h (135–136)	148	78/80	+	N.T. ⁱ
	4 (137–138)	116	69/80	+	N.T.
	1 + 2 + 3	—	80/80	+	N.T.
	1 (74–75)	167	73/80	+	N.T.
<i>P. mirabilis</i>	2 (133–134)	123	80/80	+	N.T.
	1 ^e (83–84)	139	79/80	+	N.T.
<i>P. aeruginosa</i>	2 ^e (85–86)	223	80/80	+	N.T.
	1 (98–99)	126	79/80	+	+
<i>S. saprophyticus</i>	2 (139–140)	190	80/80	+	N.T.
	1 (112–113)	157	79/80	+	N.T.
<i>M. catarrhalis</i>	2 (118–119)	118	80/80	+	N.T.
	3 (160–119)	137	80/80	+	N.T.
	1 ^e (154–155)	217	80/80	+	N.T.
<i>H. influenzae</i>	1 ^e (156–157)	134	80/80	+	N.T.
<i>S. pneumoniae</i>	2 ^e (158–159)	197	74/80	+	N.T.
	3 (78–79)	175	67/80	+	N.T.
	1 (147–148)	175	80/80	+	N.T.
<i>S. epidermidis</i>	2 (145–146)	125	80/80	+	N.T.
	1 (152–153)	108	80/80	+	N.T.
<i>S. aureus</i>	2 (149–150)	151	80/80	+	N.T.
	3 (149–151)	176	80/80	+	N.T.
	1 ^e (141–142)	213	80/80	+	N.T.
<i>S. pyogenes</i> ^f	2 ^e (143–144)	157	24/24	+	N.T.
	1 ^e (126–127)	241	194/195 ^g	+	+

See notes on next page

- [0090] a All primer pairs are specific in PCR assays since no amplification was observed with DNA from 66 different species of both Gram positive and Gram negative bacteria other than the species of interest (Table 5).
- [0091] b The ubiquity was normally tested on 80 strains of the species of interest. All retained primer pairs amplified at least 90% of the isolates. When combinations of primers were used, an ubiquity of 100% was reached.
- [0092] c For all primer pairs and multiplex combinations, PCR amplifications directly performed from a bacterial colony were 100% species-specific.
- [0093] d PCR assays performed directly from urine specimens.
- [0094] e Primer pairs derived from data bank sequences. Primer pairs with no "e" are derived from our species-specific fragments.
- [0095] f For *S. pyogenes*, primer pair #1 is specific for Group A Streptococci (GAS). Primer pair #2 is specific for the GAS-producing exotoxin A gene (SpeA).
- [0096] g Ubiquity tested on 195 isolates from 23 species representative of bacterial pathogens commonly encountered in clinical specimens.
- [0097] h Optimizations are in progress to eliminate non-specific amplification observed with some bacterial species other than the target species.
- [0098] N.T.: not tested.

TABLE 8

Selected antibiotic resistance genes for diagnostic purposes.			
Genes	Antibiotics	Bacteria ^a	SEQ ID NO
(bla _{TEM}) TEM-1	β-lactams	Enterobacteriaceae, Pseudomonadaceae, Haemophilus, Neisseria	161
(bla _{OB}) ROB-1	β-lactams	Haemophilus, Pasteurella	162
(bla _{SHV}) SHV-1	β-lactams	Klebsiella and other	163
aadB, aacC1, aacC2, aacC3, aacA4	Aminoglycosides	Enterobacteriaceae, Pseudomonadaceae	164, 165, 166 167, 168
mecA	β-lactams	Staphylococci	169
vanH, vanA, vanX	Vancomycin	Enterococci	170
satA	Macrolides	Enterococci	173
aacA-aphD	Aminoglycosides	Enterococci, Staphylococci	174
vat	Macrolides	Staphylococci	175
vga	Macrolides	Staphylococci	176
msrA	Erythromycin	Staphylococci	177
Int and Sul conserved sequences	β-lactams, trimethoprim, aminoglycosides, antiseptic, chloramphenicol	Enterobacteriaceae, Pseudomonadaceae	171, 172

^aBacteria having high incidence for the specified antibiotic resistance genes. The presence in other bacteria is not excluded.

[0099]

Annex I: Specific and ubiquitous oligonucleotides probes for hybridization

SEQ ID NO	Nucleotide Sequence	Originating DNA fragment	
		SEQ ID NO	Nucleotide position
<u>Bacterial species: <i>Escherichia coli</i></u>			
44	5'-CAC CCG CTT GCG TGG CAA GCT GCC C	5 ^a	213-237
45	5'-CGT TTG TGG ATT CCA GTT CCA TCC G	5 ^a	489-513
48	5'-TGA AGC ACT GGC CGA AAT GCT GCG T	6 ^a	759-783
49	5'-GAT GTA CAG GAT TCG TTG AAG GCT T	6 ^a	898-922
50	5'-TAG CGA AGG CGT AGC AGA AAC TAA C	7 ^a	1264-1288
51	5'-GCA ACC CGA ACT CAA CGC CGG ATT T	7 ^a	1227-1251
52	5'-ATA CAC AAG GGT CGC ATC TGC GGC C	7 ^a	1313-1337

-continued

 Annex I: Specific and ubiquitous oligonucleotides probes
 for hybridization

SEQ ID NO	Nucleotide Sequence	Originating DNA fragment	
		SEQ ID NO	Nucleotide position
53	5'-TGC GTA TGC ATT GCA GAC CTT GTG GC	7 ^a	111-136
54	5'-GCT TTC ACT GGA TAT CGC GCT TGG G	7 ^a	373-397
<u>Bacterial species: <i>Proteus mirabilis</i></u>			
70 ^b	5'-TGG TTC ACT GAC TTT GCG ATG TTT C	12	23-47
71	5'-TCG AGG ATG GCA TGC ACT AGA AAA T	12	53-77
72 ^b	5'-CGC TGA TTA GGT TTC GCT AAA ATC TTA TTA	12	80-109
73	5'-TTG ATC CTC ATT TTA TTA ATC ACA TGA CCA	12	174-203

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the Sequence listing**[0100]**

 Annex I: Specific and ubiquitous oligonucleotides probes
 for hybridization

SEQ ID NO	Nucleotide Sequence	Originating DNA fragment	
		SEQ ID NO	Nucleotide position
<u>Bacterial species: <i>Proteus mirabilis</i></u>			
76	5'-CCG CCT TTA GCA TTA ATT GGT GTT TAT AGT	13	246-275
77	5'-CCT ATT GCA GAT ACC TTA AAT GTC TTG GGC	13	291-320
80 ^b	5'-TTG AGT GAT GAT TTC ACT GAC TCC C	14	18-42
81	5'-GTG AGA CAG TGA TGG TGA GGA CAC A	15 ^a	1185-1209
82	5'-TGG TTG TCA TGC TGT TTG TGT GAA AAT	15 ^a	1224-1230
<u>Bacterial species: <i>Klebsiella pneumoniae</i></u>			
57	5'-GTG GTG TCG TTC AGG GGT TTC AC	8	45-67
58	5'-GCG ATA TTC ACA CCC TAC GCA GCC A	9	161-185
59 ^b	5'-GTC GAA AAT GCC GGA AGA GGT ATA CG	9	203-228
60 ^b	5'-ACT GAG CTG CAG ACC GGT AAA ACT CA	9	233-258
63 ^b	5'-CGT GAT GGA TAT TCT TAA CGA AGG GC	10	250-275
64 ^b	5'-ACC AAA CTG TTG AGC CGC CTG GA	10	201-223
65	5'-GTG ATC GCC CCT CAT CTG CTA CT	10	77-99
66	5'-CGC CCT TCG TTA AGA ATA TCC ATC AC	10	249-274
69	5'-CAG GAA GAT GCT GCA CCG GTT GTT G	11 ^a	296-320

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the Sequence listing

[0101]

Annex I: Specific and ubiquitous oligonucleotides probes for hybridization			
		Originating DNA fragment	
SEQ ID NO	Nucleotide Sequence	SEQ ID NO	Nucleotide position
<u>Bacterial species: <i>Pseudomonas aeruginosa</i></u>			
87	5'-AAT GCG GCT GTA CCT CGG CGC TGG T 18 ^a		2985-3009
88	5'-GGC GGA GGG CCA GTT GCA CCT GCC A 18 ^a		2929-2953
89	5'-AGC CCT GCT CCT CGG CAG CCT CTG C 18 ^a		2821-2845
90	5'-TGG CTT TTG CAA CCG CGT TCA GGT T 18 ^a		1079-1103
91	5'-GCG CCC GCG AGG GCA TGC TTC GAT G 19 ^a		705-729
92	5'-ACC TGG GCG CCA ACT ACA AGT TCT A 19 ^a		668-692
93	5'-GGC TAC GCT GCC GGG CTG CAG GCC G 19 ^a		505-529
94	5'-CCG ATC TAG ACC ATC GAG ATG GGC G 20 ^a		1211-1235
95	5'-GAG CGC GGC TAT GTG TTC GTC GGC T 20 ^a		2111-2135
<u>Bacterial species: <i>Streptococcus pneumoniae</i></u>			
120	5'-TCT GTG CTA GAG ACT GCC CCA TTT C 30		423-447
121	5'-CGA TGT CTT GAT TGA GCA GGG TTA T 31 ^a		1198-1222

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the Sequence listing

[0102]

Annex I: Specific and ubiquitous oligonucleotides probes for hybridization			
		Originating DNA fragment	
SEQ ID NO	Nucleotide Sequence	SEQ ID NO	Nucleotide position
<u>Bacterial species: <i>Staphylococcus saprophyticus</i></u>			
96	5'-CGT TTT TAC CCT TAC CTT TTC GTA CTA CC	21	45-73
97 ^b	5'-TCA GGC AGA GGT AGT ACG AAA AGG TAA GGG	21	53-82
100	5'-CAC CAA GTT TGA CAC GTG AAG ATT CAT	22	89-115
101 ^b	5'-ATG AGT GAA GCG GAG TCA GAT TAT GTG CAG	23	105-134
102	5'-CGC TCA TTA CGT ACA GTG ACA ATC G	24	20-44
103	5'-CTG GTT AGC TTG ACT CTT AAC AAT CTT GTC	24	61-90
104 ^b	5'-GAC GCG ATT GTC ACT GTA CGT AAT GAG CGA	24	19-48
<u>Bacterial species: <i>Moraxella catarrhalis</i></u>			
108	5'-GCC CCA AAA CAA TGA AAC ATA TGG T-3'	28	81-105
109	5'-CTG CAG ATT TTG GAA TCA TAT CGC C-3'	28	126-130
110	5'-TGG TTT GAC CAG TAT TTA ACG CCA T-3'	28	165-189

-continued

 Annex I: Specific and ubiquitous oligonucleotides probes
for hybridization

SEQ ID NO	Nucleotide Sequence	Originating DNA fragment	
		SEQ ID NO	Nucleotide position
111	5'-CAA CGG CAC CTG ATG TAC CTT GTA C-3'	28	232-256
114	5'-TTA CAA CCT GCA CCA CAA GTC ATC A-3'	29	97-121
115	5'-GTA CAA ACA AGC CGT CAG CGA CTT A-3'	29	139-163
116	5'-CAA TCT GCG TGT GTG CGT TCA CT-3'	29	178-200
117	5'-GCT ACT TTG TCA GCT TTA GCC ATT CA-3'	29	287-312

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the Sequence listing**[0103]**

 Annex I:
Specific and ubiquitous oligonucleotides probes
for hybridization

SEQ ID NO	Nucleotide sequence	Originating DNA fragment	
		SEQ ID NO	Nucleotide position
Bacterial species: <i>Haemophilus influenzae</i>			
105 ^b	5'-GCG TCA GAA AAA GTA GGC GAA ATG AAA G	25	138-165
106 ^b	5'-AGC GGC TCT ATC TTG TAA TGA CAC A	26 ^a	770-794
107 ^b	5'-GAA ACG TGA ACT CCC CTC TAT ATA A	27 ^a	5184-5208
Universal probes^c			
122 ^b	5'-ATC CCA CCT TAG GCG GCT GGC TCC A	-	-
123	5'-ACG TCA AGT CAT CAT GGC CCT TAC GAG TAG G	-	-
124 ^b	5'-GTG TGA CGG GCG GTG TGT ACA AGG C	-	-
125 ^b	5'-GAG TTG CAG ACT CCA ATC CGG ACT ACG A	-	-
128 ^b	5'-CCC TAT ACA TCA CCT TGC GGT TTA GCA GAG AG	-	-
129	5'-GGG GGG ACC ATC CTC CAA GGC TAA ATA C	-	-
130 ^b	5'-CGT CCA CTT TCG TGT TTG CAG AGT GCT GTG TT	-	-

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the Sequence listing^cUniversal probes were derived from 16S or 23S ribosomal RNA gene sequences not included in the Sequence listing

[0104]

<u>Annex II: Specific and ubiquitous primers for DNA amplification</u>					
			<u>Originating DNA fragment</u>		
SEQ ID NO	Nucleotide Sequence		SEQ ID NO	Nucleotide position	
<u>Bacterial species: <i>Escherichia coli</i></u>					
42	5'-GCT TTC CAG CGT CAT ATT G	4		177-195	
43 ^b	5'-GAT CTC GAC AAA ATG GTG A	4		260-278	
46	5'-TCA CCC GCT TGC GTG GC	5 ^a		212-228	
47 ^b	5'-GGA ACT GGA ATC CAC AAA C	5 ^a		490-508	
55	5'-GCA ACC CGA ACT CAA CGC C	7 ^a		1227-1245	
56 ^b	5'-GCA GAT GCG ACC CTT GTG T	7 ^a		1315-1333	
131	5'-CAG GAG TAC GGT GAT TTT TA	3		60-79	
132 ^b	5'-ATT TCT GGT TTG GTC ATA CA	3		174-193	
<u>Bacterial species: <i>Enterococcus faecalis</i></u>					
38	5'-GCA ATA CAG GGA AAA ATG TC	1 ^a		69-88	
39 ^b	5'-CTT CAT CAA ACA ATT AAC TC	1 ^a		249-268	
40	5'-GAA CAG AAG AAG CCA AAA AA	2 ^a		569-588	
41 ^b	5'-GCA ATC CCA AAT AAT ACG GT	2 ^a		670-689	

^aSequences from data banks^b These sequences are from the opposite DNA strand of the sequences given in the Sequence listing

[0105]

<u>Annex II: Specific and ubiquitous primers for DNA amplification</u>					
			<u>Originating DNA fragment</u>		
SEQ ID NO	Nucleotide Sequence		SEQ ID NO	Nucleotide position	
<u>Bacterial species: <i>Klebsiella pneumoniae</i></u>					
61	5'-GAC AGT CAG TTC GTC AGC C	9		37-55	
62 ^b	5'-CGT AGG GTG TGA ATA TCG C	9		161-179	
67	5'-TCG CCC CTC ATC TGC TAC T	10		81-99	
68 ^b	5'-GAT CGT GAT GGA TAT TCT T	10		260-278	
135	5'-GCA GCG TGG TGT CGT TCA	8		40-57	
136 ^b	5'-AGC TGG CAA CGG CTG GTC	8		170-187	
137	5'-ATT CAC ACC CTA CGC AGC CA	9		166-185	
138 ^b	5'-ATC CGG CAG CAT CTC TTT GT	9		262-281	
<u>Bacterial species: <i>Proteus mirabilis</i></u>					
74	5'-GAA ACA TCG CAA AGT CAG T	12		23-41	
75 ^b	5'-ATA AAA TGA GGA TCA AGT TC	12		170-189	

-continued

Annex II: Specific and ubiquitous primers for DNA amplification

SEQ ID NO	Nucleotide Sequence	<u>Originating DNA fragment</u>
SEQ ID NO	Nucleotide Sequence	SEQ ID NO Nucleotide position
133	5'-CGG GAG TCA GTG AAA TCA TC 14	17-36
134 ^b	5'-CTA AAA TCG CCA CAC CTC TT 14	120-139

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the Sequence listing**[0106]**

Annex II: Specific and ubiquitous primers for DNA amplification

SEQ ID NO	Nucleotide Sequence	<u>Originating DNA fragment</u>
SEQ ID NO	Nucleotide Sequence	SEQ ID NO Nucleotide position
<u>Bacterial species: <i>Staphylococcus saprophyticus</i></u>		
98	5'-CGT TTT TAC CCT TAC CTT TTC GTA CT 21	45-70
99 ^b	5'-ATC GAT CAT CAC ATT CCA TTT GTT TTT A 21	143-170
139	5'-CTG GTT AGC TTG ACT CTT AAC AAT C 24	61-85
140 ^b	5'-TCT TAA CGA TAG AAT GGA GCA ACT G 24	226-250
<u>Bacterial species: <i>Pseudomonas aeruginosa</i></u>		
83	5'-CGA GCG GGT GGT GTT CAT C 16 ^a	554-572
84 ^b	5'-CAA GTC GTG GTG GGA GGG A 16 ^a	674-692
85	5'-TCG CTG TTC ATC AAG ACC C 17 ^a	1423-1441
86 ^b	5'-CCG AGA ACC AGA CTT CAT C 17 ^a	1627-1645
<u>Bacterial species: <i>Moraxella catarrhalis</i></u>		
112	5'-GGC ACC TGA TGT ACC TTG 28	235-252
113 ^b	5'-AAC AGC TCA CAC GCA TT 28	375-391
118	5'-TGT TTT GAG CTT TTT ATT TTT TGA 29	41-64
119 ^b	5'-CGC TGA CGG CTT GTT TGT ACC A 29	137-158
160	5'-GCT CAA ATC AGG GTC AGC 29	22-39
119 ^b	5'-CGC TGA CGG CTT GTT TGT ACG A 29	137-158

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the Sequence listing

[0107]

<u>Annex II: Specific and ubiquitous primers for DNA amplification</u>					
				<u>Originating DNA fragment</u>	
SEQ ID NO	Nucleotide Sequence		SEQ ID NO	Nucleotide position	
<u>Bacterial species: Staphylococcus epidermidis</u>					
145	5'-ATC AAA AAG TTG GCG AAC CTT TTC A	36		21-45	
146 ^b	5'-CAA AAG AGC GTG GAG AAA AGT ATC A	36		121-145	
147	5'-TCT CTT TTA ATT TCA TCT TCA ATT CCA TAG	36		448-477	
148 ^b	5'-AAA CAC AAT TAC AGT CTG GTT ATC CAT ATC	36		593-622	
<u>Bacterial species: Staphylococcus aureus</u>					
149 ^b	5'-CTT CAT TTT ACG GTG ACT TCT TAG AAG ATT	37		409-438	
150	5'-TCA ACT GTA GCT TCT TTA TCC ATA CGT TGA	37		288-317	
149 ^b	5'-CTT CAT TTT ACG GTG ACT TCT TAG AAG ATT	37		409-438	
151	5'-ATA TTT TAG CTT TTC AGT TTC TAT ATC AAC	37		263-292	
152	5'-AAT CTT TGT CGG TAC ACG ATA TTC TTC ACG	37		5-34	
153 ^b	5'-CGT AAT GAG ATT TCA GTA GAT AAT ACA ACA	37		83-112	

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the sequence listing

[0108]

<u>Annex II: Specific and ubiquitous primers for DNA amplification</u>					
				<u>Originating DNA fragment</u>	
SEQ ID NO	Nucleotide Sequence		SEQ ID NO	Nucleotide position	
<u>Bacterial species: Haemophilus influenzae</u>					
154	5'-TTT AAC GAT CCT TTT ACT CCT TTT G	27 ^a		5074-5098	
155 ^b	5'-ACT GCT GTT GTA AAG AGG TTA AAA T	27 ^a		5266-5290	
<u>Bacterial species: Streptococcus pneumoniae</u>					
78	5'-AGT AAA ATG AAA TAA GAA CAG GAC AG	34		164-189	
79 ^b	5'-AAA ACA GGA TAG GAG AAC GGG AAA A	34		314-338	
156	5'-ATT TGG TGA CGG GTG ACT TT	31 ^a		1401-1420	
157 ^b	5'-GCT GAG GAT TTG TTC TTC TT	31 ^a		1515-1534	
158	5'-GAG CGG TTT CTA TGA TTG TA	35 ^a		1342-1361	
159 ^b	5'-ATC TTT CCT TTC TTG TTC TT	35 ^a		1519-1538	
<u>Bacterial species: Streptococcus pyogenes</u>					
141	5'-TGA AAA TTC TTG TAA CAG GC	32 ^a		286-305	
142 ^b	5'-GGC CAC CAG CTT GCC CAA TA	32 ^a		479-498	

-continued

Annex II: Specific and ubiquitous primers for DNA amplification

		<u>Originating DNA fragment</u>	
SEQ ID NO	Nucleotide Sequence	SEQ ID NO	Nucleotide position
143	5'-ATA TTT TCT TTA TGA GGG TG	33 ^a	966-985
144 ^b	5'-ATC CTT AAA TAA AGT TGC CA	33 ^a	1103-1122

^aSequences from data banks^bThese sequences are from the opposite DNA strand of the sequences given in the Sequence listing

[0109]

Annex II: Specific and ubiquitous primers for DNA amplification

		<u>Originating DNA fragment</u>	
SEQ ID NO	Nucleotide Sequence	SEQ ID NO	Nucleotide position
<u>Universal primers^c</u>			
126	5'-GGA GGA AGG TGG GGA TGA CG	-	-
127 ^b	5'-ATG GTG TGA CGG GCG GTG TG	-	-

^a Sequences from data banks^b These sequences are from the opposite DNA strand of the sequences given in the Sequence listing^c Universal primers were derived from the 16S ribosomal RNA gene sequence not included in the Sequence listing

[0110]

Annex III
Selection of Universal Probes by Alignment of the
Sequences of Bacterial 16S and 23S Ribosomal RNA Genes.

Reverse strand of SEQ ID NO: 122 **TGGAGCC AGCCGCCTA GGTGGGAT**

<i>Streptococcus salivarius</i>	1461	1510	TGAGGTAACC TTTGGAGCC AGCCGCCTAA GGTGGGATAG ATGANNGGGG
<i>Proteus vulgaris</i>	TAGCTTAACC	TTCGGGAGGG	CGCTTACCAC TTTGTGATTC ATGACTGGGG
<i>Pseudomonas aeruginosa</i>	TAGTCTAACC	GCAAGGGGGA	CGGTTACCAC GGAGTGATTC ATGACTGGGG
<i>Neisseria gonorrhoeae</i>	TAGGGTAACC	GCAAGGAGTC	CGCTTACCAC GGTATGCTTC ATGACTGGGG
<i>Streptococcus lactis</i>	TTGCCTAACC	GCAAGGAGGG	CGCTTCCTAA GGTAAAGCCG ATGACNNGGG

[0111]

Annex III. Selection of universal probes by alignment of the
sequences of bacterial 16S and 23S ribosomal RNA genes.

SEQ ID NO: 123	ACGTCAAGTC ATCATGGC CCTTACGAGT AGG
<i>Haemophilus influenzae</i>	1251 GGTNGGGATG ACGTCAAGTC ..ATCATGGC CCTTACGAGT AGGGCTACAC
	1300

-continued

Annex III. Selection of universal probes by alignment of the
sequences of bacterial 16S and 23S ribosomal RNA genes.

<i>Neisseria gonorrhoeae</i>	GGTGGGGATG ACGTCAAGTC ..CTCATGGC CCTTATGACC AGGGCTTCAC
<i>Pseudomonas cepacia</i>	GGTNGGGATG ACGTCAAGTC ..CTCATGGC CCTTATGGGT AGGGCTTCAC
<i>Serratia marcescens</i>	GGTGGGGATG ACGTCAAGTC ..ATCATGGC CCTTACGAGT AGGGCTACAC
<i>Escherichia coli</i>	GGTGGGGATG ACGTCAAGTC ..ATCATGGC CCTTACGACC AGGGCTACAC
<i>Proteus vulgaris</i>	GGTGGGGATG ACGTTAAGTC GTATCATGGC CCTTACGAGT AGGGCTACAC
<i>Pseudomonas aeruginosa</i>	GGTGGGGATG ACGTCAAGTC ..ATCATGGC CCTTACGGCN AGGGCTACAC
<i>Clostridium perfringens</i>	GGTGGGGATG ACGTNNAATC ..ATCATGCC CNTTATGTGT AGGGCTACAC
<i>Mycoplasma hominis</i>	GGTGGGGATG ACGTCAAATC ..ATCATGCC TCTTACGAGT GGGCCACAC
<i>Helicobacter pylori</i>	GGTGGGGACG ACGTCAAGTC ..ATCATGGC CCTTACGCCT AGGGCTACAC
<i>Mycoplasma pneumoniae</i>	GGAAGGGATG ACGTCAAATC ..ATCATGCC CCTTATGTCT AGGGCTGCAA

[0112]

Annex III. Selection of universal probes by alignment of the
sequences of bacterial 16S and 23S ribosomal RNA genes.

Reverse of the probe SEQ ID NO: 124	GCCTTGTACA CACCGCCCGT CACAC
	1451 1490
<i>Escherichia coli</i>	ACGTTCCCGG GCCTTGTACA CACCGCCCGT CACACCATGG
<i>Neisseria gonorrhoeae</i>	ACGTTCCCGG NNCTTGTACA CACCGCCCGT CACACCATGG
<i>Pseudomonas cepacia</i>	ACGTTCCCGG GTCTTGTACA CACNGCCCGT CACACCATGG
<i>Serratia marcescens</i>	ACGTTCCCGG GCCTTGTACA CACCGCCCGT CACACCATGG
<i>Proteus vulgaris</i>	ACGTTCCCGG GCCTTGTACA CACCGCCCGT CACACCATGG
<i>Haemophilus influenzae</i>	ACGTTCCCGG GCNTTGTACA CACCGCCCGT CACACCATGG
<i>Pseudomonas aeruginosa</i>	ACGTTCCCGG GCCTTGTACA CACCGCCCGT CACACCATGG
<i>Clostridium perfringens</i>	ACGTTCCCGG GTCTTGTACA CACCGCNCGT CACACCATGA
<i>Mycoplasma hominis</i>	ACGTTCTCGG GTCTTGTACA CACCGCCCGT CACACCATGG
<i>Helicobacter pylori</i>	ACGTTCCCGG GTCTTGTACT CACCGCCCGT CACACCATGG
<i>Mycoplasma pneumoniae</i>	ACGTTCTCGG GTCTTGTACA CACCGCCCGT CAAACTATGA

[0113]

Annex III. Selection of universal probes by alignment of the
sequences of bacterial 16S and 23S ribosomal RNA genes.

Reverse strand of SEQ ID NO 125:	TCG TAGTCCGGAT TGGAGTCTGC AACTC
	1361 1400
<i>Escherichia coli</i>	AAGTGCCTCG TAGTCCGGAT TGGAGTCTGC AACTCGACTC
<i>Neisseria gonorrhoeae</i>	AAACCGATCG TAGTCCGGAT TGCAGTCTGC AACTCGAGTG

-continued

 Annex III. Selection of universal probes by alignment of the sequences of bacterial 16S and 23S ribosomal RNA genes.

<i>Pseudomonas cepacia</i>	AAACCGATCG TAGTCCGGAT TGCACCTCTGC AACTCGAGTG
<i>Serratia marcescens</i>	AAGTATGTCG TAGTCCGGAT TGGAGTCTGC AACTCGACTC
<i>Proteus vulgaris</i>	AAGTCTGTCG TAGTCCGGAT TGGAGTCTGC AACTCGACTC
<i>Haemophilus influenzae</i>	AAGTACGTCT AAGTCCGGAT TGGAGTCTGC AACTCGACTC
<i>Pseudomonas aeruginosa</i>	AAACCGATCG TAGTCCGGAT CGCAGTCTGC AACTCGACTC
<i>Clostridium perfringens</i>	AAACCAGTCT CAGTTCGGAT TGTAGGCTGA AACTCGCCTA
<i>Mycoplasma hominis</i>	AAGCCGATCT CAGTTCGGAT TGGAGTCTGC AATTCGACTC
<i>Helicobacter pylori</i>	ACACC..TCT CAGTTCGGAT TGTAGGCTGC AACTCGCCTG
<i>Mycoplasma pneumoniae</i>	AAGTTGGTCT CAGTTCGGAT TGAGGGCTGC AATTCGTCCT

[0114]

 Annex III. Selection of universal probes by alignment of the sequences of bacterial 16S and 23S ribosomal RNA genes.

Reverse strand of SEQ ID NO: 128	CT CTCTGCTAAA CCGCAAGGTG ATGTATAGGG
	1991 2040
<i>Lactobacillus lactis</i>	AAACACAGCT CTCTGCTAAA CCGCAAGGTG ATGTATAGGG GTGACGCCT
<i>Escherichia coli</i>	AAACACAGCA CTGTGCAAAC ACGAAAGTGG ACGTATACGG TGTGACGCCT
<i>Pseudomonas aeruginosa</i>	AAACACAGCA CTCTGCAAAC ACGAAAGTGG ACGTATAGGG TGTGACGCCT
<i>Pseudomonas cepacia</i>	AAACACAGCA CTCTGCAAAC ACGAAAGTGG ACGTATAGGG TGTGACGCCT
<i>Bacillus stearothermophilus</i>	AAACACAGGT CTCTGCGAAG TCGTAAGGCG ACGTATAGGG GCTGACACCT
<i>Micrococcus luteus</i>	AAACACAGGT CCATGCGAAG TCGTAAGACG ATGTATATGG ACTGACTCCT
SEQ ID NO: 129	GGGGGGACC ATCCTCCAAG GCTAAATAC
	481 530
<i>Escherichia coli</i>	TGTCTGAATA TGGGGGGACC ATCCTCCAAG GCTAAATACT CCTGACTGAC
<i>Pseudomonas aeruginosa</i>	TGTCTGAACA TGGGGGGACC ATCCTCCAAG GCTAAATACT ACTGACTGAC
<i>Pseudomonas cepacia</i>	TGTCTGAAGA TGGGGGGACC ATCCTCCAAG GCTAAATACT CGTGATCGAC
<i>Lactobacillus lactis</i>	AGTTTGAATC CGGGAGGACC ATCTCCCAAC CCTAAATACT CCTTAGTGAC
<i>Micrococcus luteus</i>	CGTGTGAATC TGCCAGGACC ACCTGGTAAG CCTGAATACT ACCTGTTGAC

[0115]

 Annex III. Selection of universal probes by alignment of the sequences of bacterial 16S and 23S ribosomal RNA genes.

Reverse strand of SEQ ID NO: 130	AACACAGCA CTCTGCAAAC ACGAAAGTGG ACG
	1981 2030
<i>Pseudomonas aeruginosa</i>	TGTTTATTAA AACACAGCA CTCTGCAAAC ACGAAAGTGG ACGTATAGGG

-continued

Annex III. Selection of universal probes by alignment of the sequences of bacterial 16S and 23S ribosomal RNA genes.

<i>Escherichia coli</i>	TGTTTATTAA AAACACAGCA CTGTGCAAAC ACGAAAGTGG ACGTATACGG
<i>Pseudomonas cepacia</i>	TGTTTAATAA AAACACAGCA CTCTGCAAAC ACGAAAGTGG ACGTATAGGG
<i>Bacillus stearothermophilus</i>	TGTTTATCAA AAACACAGGT CTCTGCGAAG TCGTAAGGCG ACGTATAGGG
<i>Lactobacillus lactis</i>	TGTTTATCAA AAACACAGCT CTCTGCTAAA CCACAAGGTG ATGTATAGGG
<i>Micrococcus luteus</i>	TGTTTATCAA AAACACAGGT CCATGCGAAG TCGTAAGACG ATGTATATGG

[0116]

Annex IV. Selection of the universal PCR primers by alignment of the bacterial 16S ribosomal RNA gene

SEQ ID NO: 126	GGAGGAA GGTGGGGATG ACG		
Reverse strand of SEQ ID NO: 127		CA CACCGCCCGT CACACCAT	
	1241	1270.....1461	1490
<i>Escherichia coli</i>	ACTGGAGGAA GGTGGGGATG ACGTCAAGTC.....GCCTTGTACA	CACCGCCCGT	CACACCATGG
<i>Neisseria gonorrhoeae</i>	GCCGGAGGAA GGTGGGGATG ACGTCAAGTC.....NNCTTGTACA	CACCGCCCGT	CACACCATGG
<i>Pseudomonas cepacia</i>	ACCGGAGGAA GGTNGGGATG ACGTCAAGTC.....GTCTTGTACA	CACNGCCCGT	CACACCATGG
<i>Serratia marcescens</i>	ACTGGAGGAA GGTGGGGATG ACGTCAAGTC.....GCCTTGTACA	CACCGCCCGT	CACACCATGG
<i>Proteus vulgaris</i>	ACCGGAGGAA GGTGGGGATG ACGTCAAGTC.....GCCTTGTACA	CACCGCCCGT	CACACCATGG
<i>Haemophilus influenzae</i>	ACTGGAGGAA GGTNGGGATG ACGTCAAGTC.....GCNTTGTACA	CACCGCCCGT	CACACCATGG
<i>Legionella pneumophila</i>	ACCGGAGGAA GCGGGGATG ACGTCAAGTC.....GCCTTGTACA	CACCGCCCGT	CACACCATGG
<i>Pseudomonas aeruginos</i>	ACCGGAGGAA GGTGGGGATG ACGTCAAGTC.....GCCTTGTACA	CACCGCCCGT	CACACCATGG
<i>Clostridium perfringens</i>	CCAGGAGGAA GGTGGGGATG ACGTNNAATC.....GTCTTGTACA	CACCGCNCGT	CACACCATGA
<i>Mycoplasma hominis</i>	CTGGGAGGAA GGTGGGGATG ACGTCAAATC.....GTCTTGTACA	CACCGCCCGT	CACACCATGG
<i>Helicobacter pylori</i>	GGAGGAGGAA GGTGGGGACG ACGTCAAGTC.....GTCTTGTACT	CACCGCCCGT	CACACCATGG
<i>Mycoplasma pneumoniae</i>	ATTGGAGGAA GGAAGGGATG ACGTCAAATC.....GTCTTGTACA	CACCGCCCGT	CAAACATATGA

[0117]

SEQUENCE LISTING

<160> NUMBER OF SEQ ID NOS: 177

<210> SEQ ID NO 1

<211> LENGTH: 1817

<212> TYPE: DNA

<213> ORGANISM: Enterococcus faecalis

<400> SEQUENCE: 1

acagtaaaaa agttgttaac gaatgaattt gtaacaact tttttgctat ggtattgagt 60

tatgaggggc aatacagga aaaatgtcgg ctgattaagg aatttagata gtgccggtta 120

-continued

gtagttgtct ataatgaaaa tagcaacaaa tatttacgca gggaaagggg cggtcgttta	180
acgggaaaaa ttagggagga taaagcaata cttttgttgg gaaaagaaat aaaaggaaac	240
tggggaagga gtaattgtt tgatgaaggg aaataaaatt ttatacattt taggtacagg	300
catctttgtt ggaagttcat gtctattttc ttcacttttt gtgcccgcag aagaacaagt	360
ttattcagaa agtgaagttt caacagtttt atcgaagttg gaaaaggagg caatttctga	420
ggcagctgct gaacaatata cggttgtaga tcgaaaagaa gacgcgtggg ggatgaagca	480
tcttaagtta gaaaagcaaa cggaaggcgt tactgttgat tcagataatg tgattattca	540
tttagataaa aacgggtgcag taacaagtggt tacaggaaat ccagttgatc aagttgtgaa	600
aattcaatcg gttgatgcaa tcggtgaaga aggagttaaa aaaattgttg cttctgataa	660
tccagaaact aaagatcttg tcttttttagc tattgacaaa cgtgtaaata atgaagggca	720
attattttat aaagtcagag taacttcttc accaactggt gaccccgtat cattggttta	780
taaagtgaac gtcacagatg gaacaattat ggaaaaacaa gatttaacgg aacatgtcgg	840
tagtgaagta acgttaaaaa actcttttca agtaacgttt aatgtaccag ttgaaaaaag	900
caatacggga attgctttac acggaacgga taacacaggg gtttaccatg cagtagttga	960
tggcaaaaat aattattcta ttattcaagc gccatcacta gcgacattaa atcagaatgc	1020
tattgacgcc tatacgcagt gaaaatttgt gaaaacatat tatgaagatc atttccaacg	1080
acacagtatt gatgatcgag ggatgccat cttgtcagtt gttgatgaac aacatccaga	1140
tgcttatgac aatgcttttt gggatggaaa agcaatgcgt tatggtgaaa caagtacacc	1200
aacaggaaaa acgtatgctt cctcttttaga tgtagttggt catgaaatga cacatggtgt	1260
gacggaacat actgccgggt tagaatattt aggacaatca ggtgccttga atgaatctta	1320
ttctgatttg atgggttata ttatttcggg tgcattctaat ccagaaattg gtgcggatac	1380
tcagagtgtt gaccgaaaaa caggtattcg aaatttacia acgccaagta aacacggaca	1440
accagaaacc atggctcaat acgacgatcg agcacggtat aaaggaaacgc cttattatga	1500
tcaaggcggg gttcattata acagtggaat tattaatcgg attggttaca ccattatcca	1560
gaacttaggc attgaaaaag cacagactat tttctacagc tcgttagtaa attacttaac	1620
acctaaagca caattcagtg atgctcgtga tgcgatgctt gctgctgcaa aagttcaata	1680
tggcgatgaa gcagcttcag tgggtgcagc agcctttaac tctgctggaa tcggagctaa	1740
agaagacatt caggtaaacc aaccaagtga atctgttctg gtcaatgaat gaaaaaaatt	1800
ccccaattaa ataaaaa	1817

<210> SEQ ID NO 2

<211> LENGTH: 2275

<212> TYPE: DNA

<213> ORGANISM: Enterococcus faecalis

<400> SEQUENCE: 2

ggtaccaaaag aaaaaaacga acgccacaac caacagcctc taaagcaaca cctgcttctg	60
aaattgaggg agatttagca aatgtcaatg agattctttt ggttcacgat gatcgtgtcg	120
ggtcagcaac gatgggaatg aaagtcttag aagaaatttt agataaagag aaaatttcaa	180
tgccgattcg aaaaattaat attaataaat taactcaaca aacacaggct ttaattgtca	240
caaaagctga actaacggaa caagcacgta aaaaagcacc gaaagcgaca cacttatcag	300

-continued

taaaaagtta	tggttaatcc	ccaaaaatat	gaaacagtgg	gtttcgctct	taaaagaaag	360
tgcc>tagaga	ggaagaaaac	aatggaaaat	cttacgaata	tttcaattga	attaatacaa	420
cagtttaata	caaaagaaga	agctattcgc	ttttccggcc	agaaactagt	cgaggcaggc	480
tgtgttgagc	ccgcttatat	cgaagcaatg	attgaaagag	accaattgct	atctgccccat	540
atggggaatt	ttattgccat	tcctcatgga	acagaagaag	ccaaaaaatt	agtgaaaaaa	600
tcaggaatct	gtgtagtgca	agtcccagag	ggcgtaatt	ttggcaccga	agaagatgaa	660
aaaattgcta	ccgtattatt	tgggattgcc	ggagtcggtg	aagaacattt	gcaattagtc	720
caacaaattg	cactttattg	tagtgatatg	gataacgtgg	tgcaacttgc	cgatgcatta	780
agtaaagaag	aaataacaga	aaatttagcc	attgcttaaa	ggagagaata	agaatgaacg	840
cagtacattt	tggagcagga	aatattggac	gcggttttat	tggcgaaatt	ttagctaaaa	900
cgggtttcat	attaccgttt	gtggatgta	atggaaacca	tcatcaagcg	ttaaagaac	960
gtaaaagtta	tacaattgaa	ttggccgatg	cctcacatca	acaaattaac	gttgaaaatg	1020
tgaccgggtt	aaataacatg	acagaaccag	aaaaagtagt	agaagcaatt	gcggaagccg	1080
athtagtcac	gacggcaatt	ggtcctaata	ttttaccaag	aattgctgaa	ttaattgctc	1140
aagggaattg	tgcacgtgcc	gaagcaaatt	gtcaaaacgg	cccgtggat	attatcgctt	1200
gtgaaaatat	gattggtggt	tcaacctttt	tagcagaaga	agtggccata	atatttgaaa	1260
aaccagctt	atctgaacaa	tggattggtt	ttcctgatgc	ggcagttgat	cggtattgtc	1320
cattacaaaa	acataaagat	ccactttttg	ttcaagttga	gcctttttgt	gaatgggtca	1380
ttgatgatac	caaccgaaaa	gccaaagaga	ttcagttaga	aggcgtcatt	acttgtcgat	1440
tagagccgta	tattgaacga	aaattattta	gtgtaaccag	tggccatgct	acagttgcct	1500
atacaggggc	gtgttaggac	tatcaaacca	ttgacgaagc	gatgcaggac	gccttagtgg	1560
tagcgcaact	caaatcagtt	ttgcaggaaa	ccggtaaact	tttagtggcc	aaatggaatt	1620
ttgatgaaca	agaacatgca	gcctatattg	aaaaaattat	caaccgtttc	caaaataaat	1680
atatttcaga	tgctattaca	cgtgtagcac	ggacaccaat	cagaaaatta	ggtgcgcaag	1740
aacggtttat	tcgaccaatc	cgtgaattac	aggaacgcaa	tctagtgtcg	gccgcattta	1800
tagcaatgat	tggtattgtc	tttaattatc	atgatccaga	agatgaacaa	agccgtcaat	1860
tacaggaaat	gcttgaccaa	gaaagtgttg	atacagtgga	tcgctgaagt	aacgggcatt	1920
gaagatccag	aaacggttaa	aaatattaa	caaacgtag	aactgctatg	cgcgaccaca	1980
agtagcataa	ttaacaaaat	ccttctacca	agatacttca	catttcttaa	ttaaagaaaa	2040
aacaaccgcg	cctcacctga	gccgaccccc	aaaagttaga	cctagaaatc	taacttttgg	2100
aggttttttt	gtatggcaaa	atacagtttt	gaaatttaaa	cttaaacttg	ttcatgacta	2160
cttatatggt	caaggaggtc	taaggtttct	cgcaaagaag	tatgggttta	aagatagtct	2220
caataagca	aatggataaa	tgcctataaa	gaacttggtg	aagaaggggg	gatcc	2275

<210> SEQ ID NO 3

<211> LENGTH: 227

<212> TYPE: DNA

<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 3

gatccgccat ggggtgtttt ccgattgagg attttataga tggtttctgg cgacctgcac 60

-continued

```
aggagtacgg tgatttttaa ttattgcaat tgcacaagag tcagttctcc cccaaagaca 120
gcaccggtat caatataatg caggttgcca atatccacgc gatggcgcaa aggtgtatga 180
ccaaaccaga aatgatcggc cacctgcacg gccagttcgc gagtcgg 227
```

```
<210> SEQ ID NO 4
<211> LENGTH: 278
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli
```

```
<400> SEQUENCE: 4
```

```
gatctaaatc aaattaattg gttaaagata accacagcgg gcccgacata aactctgaca 60
agaagttaac aacctataaa cctgcacagg acgcgaacat gtcttctcat ccgtatgtca 120
cccagcaaaa taccgccgctg gcggaacgaca ccaactctgat gtccactacc gatctcgttt 180
tccagcgtca tattggggcg cgctacgttg gggcgtgggc gtaattggtc aatcaggcgc 240
ggggtcagcg gataaacatt caccattttg tcgagatc 278
```

```
<210> SEQ ID NO 5
<211> LENGTH: 1596
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli
```

```
<400> SEQUENCE: 5
```

```
atggctgaca ttctgctgct cgataaatc gactctttta cgtacaacct gccagatcag 60
ttgcgagca atgggcaata cgtggtgatt taccgcaacc atataccggc gcaaacctta 120
attgaacgct tggcgaccat gagtaatccg gtgctgatgc tttctcctgg ccccggtgtg 180
ccgagcgaag ccggttgat gccggaactc ctcaccgct tgcgtggcaa gctgccatt 240
attggcattt gcctcggaca tcaggcgatt gtcgaagctt acgggggcta tgtcggtcag 300
gcgggcgaaa ttctccacgg taaagcctcc agcattgaac atgacggta gccgatgttt 360
gccgattaa caaacccgct gccggtggcg cgttatcact cgctggttgg cagtaacatt 420
ccggccggtt taaccatcaa cccccattt aatggcatgg tgatggcagt acgtcacgat 480
gcggatcgcg tttgtgatt ccagttccat ccggaatcca ttctcaccac ccagggcgt 540
cgctgctgg aacaaacgct gccctggcg cagcataaac tagagccagc caacacgctg 600
caaccgattc tggaaaaact gtatcaggcg cagacgctta gccaaacaaga aagccaccag 660
ctgttttcag cgggtgtgcg tggcgagctg aagccggaac aactggcggc gccgctgtg 720
agcatgaaa ttcgcggtga gcaccogaac gagatcgcg gggcagcaac cgcgctactg 780
gaaaacgca cgccttccc gcgccggat tatctgtttg ctgatatcgt cggtaactggc 840
ggtgacggca gcaacagat caatatttct accgccagtg cgtttgtcgc cgcggcctgt 900
gggctgaaag tggcgaaca cggcaaccgt agcgtctcca gtaaatctgg ttcgtccgat 960
ctgctggcgg cgttcggtat taatcttgat atgaacgccg ataaatcgcg ccagggcgtg 1020
gatgagttag gtgtatgttt cctctttgcg cogaagtatc acaccgatt ccgccacgcg 1080
atgccggttc gccagcaact gaaaaccgc accctgttca atgtgctggg gccattgatt 1140
aaccggcgc atccgccgct gccgttaatt ggtgtttata gtcgggaact ggtgctgccg 1200
attgccgaaa ccttgccgct gctggggtat caacgcgcgg cgggtgtgca cagcggcggg 1260
atggatgaag tttcattaca cgcgcgaca atcgttgccg aactgcatga cggcgaatt 1320
```

-continued

aaaagctatc agctcaccgc agaagacttt ggcctgacac cctaccacca ggagcaactg	1380
gcagggcgaa caccggaaga aaaccgtgac attttaacac gtttgttaca aggtaaaggc	1440
gacgcccgcc atgaagcagc cgtcgtcgcc aacgtcgcca tgtaaatgcg cctgcatggc	1500
catgaagatc tgcaagccaa tgcgcaaacc gttcttgagg tactgcgagc tggttccgct	1560
tacgacagag tcaccgcact ggcggcacga gggtaa	1596

<210> SEQ ID NO 6

<211> LENGTH: 2703

<212> TYPE: DNA

<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 6

gacgacttag ttttgacgga atcagcatag ttaatcactt cactgtggaa aatgaggaaa	60
tattatTTTT tttgcgcttc gtaattaatg gttataaggt cggccagaaa cttttcta	120
gcaagcgatg acgttttttt atgtgtctga atttgcactg tgcacaatt ccaaatcttt	180
attaacaact cacctaaaac gacgctgac cagcgtgaat actggtttcc cttatgttca	240
tcagattcat ttaagcaagg gtttcttctt cattcctgat gaaagtgcc tctaaaaaga	300
tgatcttaat aaatctatta agaatgagat ggagcacact ggatatttta cttatgaaac	360
tgtttcaact ctttacttaa tttatagagt taccttccgc tttttgaaa tacgcaacgg	420
ccattttttg cacttagata cagattttct gcgctgtatt gcattgattt gatgcta	480
ctgtggtttg cactagcttt aagtggttga gatcacattt ccttgcctcat ccccgca	540
cctccctgcc taatccccgc caggatgagg aaggtaaca tcgagcctgg caaactagc	600
ataacgttgt gttgaaaatc taagaaaagt ggaactccta tgcacaacc tattttta	660
gataagcaat ttcaggaagc gctttcacgt cagtggcagc gttatggctt aaattctg	720
gctgaaatga ctctcgcca gtggtggcta gcagtgagtg aagcactggc cgaatgct	780
cgtgctcagc cattcgccaa gccggtggcg aatcagcgac atgttaacta catctca	840
gagtttttga ttggtcgcct gacgggcaac aacctgttga atctcgctg gtatcagg	900
gtacaggatt cgttgaagcg ttatgacatc aatctgacgg acctgctgga agaagag	960
gacccggcgc tgggtaacgg tggcttggga cgtctggcgg cgtgcttctc cgactca	1020
gcaactgtcg gtcagctcgc gacgggttac ggtctgaact atcaatatg tttgttcc	1080
cagtcttttg tcgatggcaa acaggttgaa gcgccgatg actggcatcg cagtaact	1140
ccgtggttcc gccacaacga agcactgat gtgcaggtag ggattggcgg taaagtga	1200
aaagacggac gctgggagcc ggagtttacc attaccggtc aagcgtggga tctcccgt	1260
gtcggctatc gtaatggcgt ggcgcagccg ctgctctgtt ggcaggcgac gcaocgc	1320
ccgtttgatc tgactaaatt taacgacggt gatttcttgc gtgccgaaca gcaggc	1380
aatgcgga aactgaccaa agttctctat ccaaacgaca accatactgc cggtaaaa	1440
ctgocctga tgcagcaata cttccagtgt gctgttccg tagcggatat tttgctc	1500
catcatctgg cggggcgtga actgcacgaa ctggcggatt actaagtat tcagctga	1560
gataccacc caactatcgc gattccagaa ctgctgctgg tgctgatoga tgagcacc	1620
atgagctgg atgacgcttg gccattacc agcaaaactt tcgcttacac caaccata	1680
ctgatgccag aagcgttgg aacgctggat gtgaaactgg tgaaggctt actgccgc	1740

-continued

```

cacatgcaga ttattaacga aattaatact cgctttaaaa cgctggtaga gaaaacctgg 1800
ccgggcgatg aaaaagtgtg ggccaaactg cgggtggtgc acgacaaaca agtgcataatg 1860
gcgaacctgt gtgtggttgg cggtttcgcg gtgaacggtg ttgcggcgct gcaactcggat 1920
ctggtggtga aagatctgtt cccggaatat caccagctat ggccgaacaa attccataac 1980
gtcaccaacg gtattacccc acgtcgtctg atcaaacagt gcaacccggc actggcggct 2040
ctgttgata aatcactgca aaaagagtgg gctaacgatc tcgatcagct gatcaatctg 2100
gttaaattgg ctgatgatgc gaaattccgt cagctttatc gcgtgatcaa gcaggcgaat 2160
aaagtccgtc tggcggagt tgtgaaagtt cgtaccggtg ttgacatcaa tccacaggcg 2220
attttcgata ttcagatcaa acgtttgac gagtacaac gccagcacct gaatctgctg 2280
cgtattctg cgttgtaca agaaattcgt gaaaaccgc aggctgatcg cgtaccgcg 2340
gtcttcctc tcggcgcgaa agcggcaccg ggctactacc tggctaagaa tattatcttt 2400
gcgatcaaca aagtggctga cgtgatcaac aacgatccgc tggttggcga taagttgaag 2460
gtggtgttcc tccgggatta ttgcgtttcg gcggcggaaa aactgatccc ggcggcggat 2520
atctccgaac aaatttcgac tgcaggtaaa gaagcttccg gtaccgcaa tatgaaactg 2580
gcgctcaatg gtgcgcttac tgcgggtacg ctggatgggg cgaacgttga aatcgcggag 2640
aaagtcggtg aagaaaatat ctttattttt ggtcatacgg tcaacaagt gaaggcaatc 2700
gac 2703

```

<210> SEQ ID NO 7

<211> LENGTH: 1391

<212> TYPE: DNA

<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 7

```

agagaagcct gtcggcaccg tctggtttgc ttttgccact gcccgcggtg aaggcattac 60
ccggcgggat gcttcagcgg cgaccgtgat cgggtgctgc gtcaggctac tgcgtatgca 120
ttgcagacct tgtggcaaca atttctaca aacacttgat actgatgag catacagtat 180
aattgcttca acagaacata ttgactatcc ggtattaccc ggcatgacag gagtaaaaat 240
ggctatcgac gaaaacaaac agaaagcgtt ggccgcagca ctgggccaga ttgagaaaca 300
atthgtgtaa ggctccatca tgcgcctggg tgaagaccgt tccatggatg tggaaacct 360
ctctaccggt tcgcttttca tggatatacg gcttggggca ggtggtctgc cgatgggccc 420
tatcgtcgaa atctacggac cggaatcttc cggtaaaacc acgctgacgc tgcaggtgat 480
cgccgcagcg cagcgtgaag gtaaaacctg tgcgtttatc gatgctgaac acgcgctgga 540
cccaatctac gcacgtaaac tgggcgtcga taticgacaac ctgctgtgct cccagccgga 600
caccgcgag caggcactgg aaatctgtga cgccctggcg cgttctggcg cagtagacgt 660
tatcgtcgtt gactcctgg cggcactgac gccgaaagcg gaaatcgaag gcgaaatcgg 720
cgactctcac atgggacctg cggcaogtat gatgagccag gcgatgcgta agctggcggg 780
taacctgaag cagtccaaca cgctgctgat cttcatcaac cagatccgta tgaaaattgg 840
tgtgatgttc ggtaaccctg aaaccactac cggtggtaac gcgctgaaat tctacgcctc 900
tgttcgtctc gacatccgtc gtatcggcgc ggtgaaagag gccgaaaacg tggtggttag 960
cgaacccgc gtgaaagtgg tgaagaacaa aatcgcctgc ccgtttaaac aggctgaatt 1020

```

-continued

```

ccagatcctc tacggcgaag gtatcaactt ctacggcgaa ctggttgacc tgggcgtaaa 1080
agagaagctg atcgagaaag caggcgcgtg gtacagctac aaaggtgaga agatcgggtca 1140
gggtaaagcg aatgcgactg cctggctgaa agataaccg gaaaccgga aagagatcga 1200
gaagaaagta cgtgagttgc tgctgagcaa cccgaactca acgccggatt tctctgtaga 1260
tgatagcgaa ggcgtagcag aaactaacga agatttttaa tcgtcttggt tgatacacia 1320
gggtcgcctc tgcggccctt ttgctttttt aagttgtaag gatatgcat gacagaatca 1380
acatcccgtc g 1391

```

```

<210> SEQ ID NO 8
<211> LENGTH: 238
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

```

```
<400> SEQUENCE: 8
```

```

tcgccaggaa ggcggcattc ggctgggtca gagtgacctg cagcgtggtg tcgttcagcg 60
ctttcacccc caacgctctg ggtccctttt gcccgagggc aatctcgcgg gcgttggcga 120
tatgcatatt gccagggtag ctccgctagg gggaggctgt tgccggcgag accagccgtt 180
gccagctcca gacgatatcc tgcgctgtaa tggccgtgcc gtcagaccag gtcagacc 238

```

```

<210> SEQ ID NO 9
<211> LENGTH: 385
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

```

```
<400> SEQUENCE: 9
```

```

cagcgtaatg cgccgcggca taacggcgcc actatcgaca gtcagttcgt cagcctgcag 60
cctgggctga atctgggacc atggcgctg ccgaactaca gcacctatag ccacagcgat 120
aacaacagcc gctgggagtc ggtttactcc tatcttgccc gcgatattca caccctacgc 180
agccagctgg tggctcgtaa tacgtatacc tcttccggca ttttcgacag tttgagtttt 240
accggtctgc agctcagttc gacaaagaga tgctgccgga tagcctgcat gctttgcgcc 300
gacgattcga gggatcgcgc gcaccaccgc ggaggctctg gtttatcaga atggttacag 360
catttataaa accaccgtcg ctacc 385

```

```

<210> SEQ ID NO 10
<211> LENGTH: 462
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

```

```
<400> SEQUENCE: 10
```

```

ctctatattc aggacgaaca tatctggacc tctggcgggg tcagttccgg ctttgatcgc 60
cctgcacccg cagcgggtga tcgcccctca tctgctactg cggcgctgca acagggcagc 120
atcgatgacg ttattcctgg ccagcaaaca gcagaccaat taaggctctga tagtggctct 180
cttctccggc cgcgcgacgg tccagggcgc tcaacagttt ggtgcatagc gctttgcggt 240
tgagatgacg cccttcgtta agaatatcca tcacgatctc cgtccatgga gagtagcgtt 300
tattccagaa tagggttttt caggatctca tggatctgcg cctgcttata gctattttgt 360
aaccagatcg cataaagtgg acgggataac gtagcgtctg ccatgaccgt atgtaacca 420
tgcttctctt tcgcccagcg agcaggtagc caacagcagc cg 462

```

-continued

<210> SEQ ID NO 11
<211> LENGTH: 730
<212> TYPE: DNA
<213> ORGANISM: *Klebsiella pneumoniae*

<400> SEQUENCE: 11

gctgaccgct aaactggggt acccgatcac tgacgatctg gacatctaca cccgtctggg	60
cggcatgggt tggcgcgctg actccaaagg caactacgct tcaaccggcg tttcccgtag	120
cgaacacgac actggcggtt ccccagtatt tgctggcggc gtagagtggg ctgttactcg	180
tgacatcgct acccgtctgg aataccagtg ggtaacaac atcggcgacg cgggcactgt	240
gggtaccgct cctgataacg gcatgctgag cctgggctt tcctaccgct tccgtcagga	300
agatgctgca ccggttggtg ctccggctcc ggctccggct ccggaagtgg ctaccaagca	360
cttaccctg aagtctgacg ttctgttcaa cttcaacaaa gctaccctga aaccggaagg	420
tcagcaggct ctggatcagc tgtacactca gctgagcaac atggatccga aagacggttc	480
cgctgttggt ctgggctaca ccgaccgat cggttccgaa gcttacaacc agcagctgtc	540
tgagaaacgt gctcagtcgg ttggtgacta cctgggtgct aaaggcatcc cgctggcaa	600
aatctccgct cgcggcatgg gtgaatccaa cccggttact ggcaacacct gtgacaacgt	660
gaaagctcgc gctgccctga tcgattgcct ggctccggat cgctcgtgtag agatcgaagt	720
taaagttatc	730

<210> SEQ ID NO 12
<211> LENGTH: 225
<212> TYPE: DNA
<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 12

cgctactggt taaatctcat ttgaaacatc gcaaagtcag tgaaccacat attcgaggat	60
ggcatgcact agaaaatatt aataagatt tagcgaacc taatcagcgc aatatcgctt	120
aattatttta ggtatgctt cttctatcct acagtcacga ggcagtgctg aacttgatcc	180
tcattttatt aatcacatga ccaatggtat aagcgtcgtc acata	225

<210> SEQ ID NO 13
<211> LENGTH: 402
<212> TYPE: DNA
<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 13

acattttaaa taggaagcca cctgataaca tccccgcagt tggatcatca gatttatagc	60
ggcatttggt atccgctaga taaaagcagt ccaacgatcc cgccaattgt tagatgaaat	120
tgactattc tttttatttg ctccgcttta tcacagtggt tttcgtttg ccgcccctgt	180
gcgccaacag ctaagaacac gcacgctctt taatgtgta ggccattaa ttaatccagc	240
gcgtccgccc tttagcatta attggtggtt atagtcctga attattaatg cctattgcag	300
atacctaaa tgtcttgggc tacaaacgtg cggcagtggt ccatagtggt ggaatggatg	360
aagtgctcatt acatgctccc acacaagtgg ctgagttaca ca	402

<210> SEQ ID NO 14
<211> LENGTH: 157
<212> TYPE: DNA

-continued

<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 14

ctgaaacgca tttatgctgg agtcagtga atcatcactc aattttcacc cgatgtattt	60
tctgttgaac aagtctttat ggcaaaaaat gcagactcag cattaaaatt aggccaagca	120
agaggtgtgg cgattttagc ggcagtcaat aatgatc	157

<210> SEQ ID NO 15

<211> LENGTH: 1348

<212> TYPE: DNA

<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 15

tttctcttta aaatcaattc ttaaagaaat tattaataat taacttgata ctgtatgatt	60
atacagtata atgagtttca acaagcaaaa tcatatacgt tttaatggta gtgacccatc	120
tttatgcttc actgcccaga gggagataac atggctattg atgaaaacaa acaaaaagca	180
ttggccgcag cacttggtca aattgaaaag caatttgta aagttctat catgcgtctg	240
ggcgaagacc gttccatgaa cgtagaaact atctctacag gatctttatc attagacgtt	300
gctttagggt cagggtgatt gccacgtggc cgtattgttg aaatctatgg ccctgaatct	360
tctgttaaaa caacctgac tctacaagt attgcctctg ctcagcgtga aggaaaaatt	420
tgtagcattta ttgatgctga acatgcatta gacccaattt atgctcaaaa gctaggtgtc	480
gatatcgata atctactctg ctctcaacct gacacagggt aacaagctct ggaatttgt	540
gatgcattat ctgcctctgg tgcggtcgat gttattgtcg tggactccgt ggcagcatta	600
acacaaaag ctgaaattga aggtgaaatt ggtgattcac acgttggtt agccgcacgt	660
atgatgagcc aagctatgag taaactagcg ggtaacctta aaaactctaa tacactgctg	720
attttcatta accaaattcg tatgaaaatc ggtgttatgt ttgtaacctt agaaccacg	780
accggtggta atgctcttaa attctatgct tctgttcggt tagacattcg tcgcattggc	840
tctgtcaaaa atgggtatga agtcattggt agtgagactc gcgttaaagt tgttaaaaat	900
aaagtggctg caccgtttaa acaagctgaa ttccaaatta tgtacgggtga aggtattaat	960
acctatggcg aactgattga tttagggtgtt aaacataagt tagtagagaa agcagggtgct	1020
tggtatagct acaatggcga aaaaattggt caaggtaaag ctaacgcaac caattactta	1080
aaagaacatc ctgaaatgta caatgagtta aacactaaat tgcgtgaaat gttgttaaat	1140
catgctggty aattcacaag tgctgcggat tttgcaggty aagagtcaga cagtgatgct	1200
gacgacacaa aagagtaatt agctggttgt catgctgttt gtgtgaaaat agaccttaa	1260
tcattggcta ttatcagcag agcatccat agaataactt gttgtataa attttattca	1320
gatggcaaaag gaagccttaa aaaagctt	1348

<210> SEQ ID NO 16

<211> LENGTH: 2167

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 16

ggtaccgctg gccgagcatc tgctcgtatc ccaccagccg ggcgacggga actgcacgat	60
ctacctggcg agcctggagc acgagcgggt tcgcttcgta cggcgtgag cgacagtcac	120

-continued

aggagaggaa acggatggga tcgcaccagg agcggccgct gatcggcctg ctgttctccg	180
aaaccggcgt caccgccgat atcgagcgct cgcacgcgta tggcgcattg ctgcgggtcg	240
agcaactgaa ccgcgagggc ggcgtcggcg gtcgcccgat cgaaacgctg tcccaggacc	300
ccggcggcga cccggaccgc tatcggctgt gcgccgagga cttcattcgc aaccgggggg	360
tacggttcct cgtgggctgc tacatgtcgc acacgcgcaa ggcggtgatg ccggtggctg	420
agcgcgccga cgcgctgctc tgctaccoga cccctacga gggcttcgag tattcgcgca	480
acatcgtcta cggcggctccg gcgccgaacc agaacagtgc gccgctggcg gcgtacctga	540
ttcgccacta cggcgcgagg gtggtgttca tcggctcggga ctacatctat ccgcgggaaa	600
gcaaccatgt gatgcgccac ctgtatcgcc agcacggcgg cacggtgctc gaggaaatct	660
acattccgct gtatccctcc gacgaogact tgcagcgcgc cgtcgcgagc atctaccagg	720
cgcgcgccga cgtggtcttc tccaccgtgg tgggcaccgg caccgccgag ctgtatcgcg	780
ccatcgcgcc tcgctacggc gacggcaggc ggccgccgat cgcagcctg accaccagcg	840
aggcggagggt gccgaagatg gagagtgcg tggcagaggg gcaggtggtg gtcgcgcctt	900
acttctccag catcgatacg cccgccagcc gggccttcgt ccaggcctgc catggtttct	960
tcccggagaa cgcgaccatc accgcctggg ccgagggcgc ctactggcag accttgttgc	1020
tcggccgcgc cgcgcaggcc gcaggcaact ggcgggtgga agacgtgcag cggcacctgt	1080
acgacatcga catcgacgcg ccacaggggc cggtcgggtt ggagcgcag aacaaccaca	1140
gccgcctgct ttcgcgcatc gcggaaatcg atgcgcgcgg cgtgttcag gtcgcgtggc	1200
agtcgccga accgattcgc cccgaccctt atgtcgtcgt gcataacctc gacgactggt	1260
ccgccagcat gggcggggga ccgctcccat gagcgcaac tcgctgctcg gcagcctgcg	1320
cgagttgcag gtgctggtcc tcaaccgcc gggggaggtc agcgcgccc tggctctgca	1380
gtgatccgc atcgggtgtt cggtcgcca gtgctggccg ccgccggaag ccttcgacgt	1440
gccggtggac gtggtcttca ccagcatttt ccagaatggc caccacgacg agatcgtgctc	1500
gtgctcgc gccgggactc cgcgcactac cctggtggcg ctggtggagt acgaaagccc	1560
cgcggtgctc tcgcagatca tcgagctgga gtgccacggc gtgatcaacc agcgcctcga	1620
tgcccaccgg gtgctgcctg tgctggtatc ggcgcggcgc atcagcgagg aaatggcgaa	1680
gtgaaagcag aagaccgagc agctccagga ccgcatcgcc ggccaggccc ggatcaacca	1740
ggccaagggt ttgctgatgc agcgcgatg ctgggacgag cgcgaggcgc accagacct	1800
gtcgcgggaa gcgatgaagc ggcgcgagcc gatcctgaag atcgcctcagg agttgctggg	1860
aaacgagccg tccgcctgag cgatccgggc cgaccagaac aataacaaga ggggtatcgt	1920
catcatgctg ggactggttc tgctgtactg tggcgcggtg ctgtttctca atgocgtctg	1980
ggtgctgggc aagatcagcg gtcgggagggt ggcggtgatc aacttcctgg tcggcgtgct	2040
gagcgcctgc gtcgcgttct acctgatctt ttccgcagca gccgggcagg gctcgtgaa	2100
ggccggagcg ctgaccctgc tattcgcttt tacctatctg tgggtggccg ccaaccagtt	2160
cctcag	2167

<210> SEQ ID NO 17

<211> LENGTH: 1872

<212> TYPE: DNA

<213> ORGANISM: Pseudomonas aeruginosa

-continued

<400> SEQUENCE: 17

```

gaattcccgg gagttcccga cgcagccacc cccaaaacac tgctaaggga gcgcctcgca    60
gggctcctga ggagatagac catgcccattt ggcaagccac tggtagggcac cttgctcgcc    120
tcgctgacgc tgctgggctt ggccaccgct cacgcccaagg acgacatgaa agccgcccggag    180
caataccagg gtgccgcttc cgcctcgat cccgctcacg tggtagcgcac caacggcgct    240
cccgacatga gtgaaagcga gttcaacgag gccaaagcaga tctacttcca acgctcgccc    300
ggttgccacg gcgtcctcgc caaggggccc accggcaagc cgctgacccc ggacatcacc    360
cagcaacgcg gccagcaata cctggaagcg ctgatcacct acggcacccc gctgggcatg    420
ccgaactggg gcagctccgg cgagctgagc aaggaacaga tcaccctgat ggccaagtac    480
atccagcaca ccccgcgcga accgcccggag tggggcatgc cggagatgcg cgaatcgtgg    540
aaggtagctgg tgaagccgga ggaccggcgg aagaaacagc tcaacgacct cgacctgccc    600
aaactgttct cggtagacct gcgcgacgcc gggcagatcg ccctggtcga cggcgacagc    660
aaaaagatcg tcaaggtcat cgataccggc tatgccgtgc atatctcgcg gatgtccgct    720
tccggccgct acctgctggt gatcggcccgc gacgcgcgga tcgacatgat cgacctgtgg    780
gccaaggagc cgaccaaggt cgcagagatc aagatcggca tcgagggcgcg ctcgggtggaa    840
agctccaagt tcaagggcta cgaggaccgc tacaccatcg ccggcgccta ctggccgccc    900
cagttcgcga tcatggacgg cgagaccctg gaaccgaagc agatcgtctc caccgcgggc    960
atgacctgag acaccagac ctaccaccgg gaaccgcgcg tggcggcgat catcgccctc    1020
cacgagcacc ccgagttcat cgtcaacgtg aaggagaccg gcaaggtcct gctgggtcaac    1080
tacaaggata tcgacaacct caccgtcacc agcatcgtg cggcgccggt cctccacgac    1140
ggcggtggg acagcagcca ccgctacttc atgaccgccc ccaacaactc caacaaggtt    1200
gccgtgatcg actccaagga ccgctgcctg tcggccctgg tcgacgtcgg caagaccccg    1260
caccgggggc gtggcgccaa cttcgtgcat cccaagtacg gcccggtgtg gagcaccagc    1320
cacctgggag acggcagcat ctcgctgacg ggcaccgac cgaagaacca tccgcagtac    1380
gcctggaaga aagtcgcccga actacagggc cagggcgggc gctcgtgtgt catcaagacc    1440
catccgaagt cctcgcacct ctacgtcgac accaccttca accccgacgc caggatcagc    1500
cagagcgtcg cgggtgtcga cctgaagaac ctcgacgcca agtaccagggt gctgccgate    1560
gccgaatggg ccgatctcgg cgaaggcgcc aagcgggtgg tgcagcccga gtacaacaag    1620
cgcggcgatg aagtctggtt ctcgggtgtg aacggcaaga acgacagctc cgcgctggtg    1680
gtggtggagc acaagacctt gaagctcaag gccgtggtca aggacccgag gctgateacc    1740
ccgaccggta agttcaacgt ctacaacacc cagcacgacg tgtactgaga cccgctgagc    1800
gggcacgccc cgcacgctcc cccctacgag gaaccgtgat gaaaccgtac gactgcttt    1860
cgtgctcgc ca                                     1872

```

<210> SEQ ID NO 18

<211> LENGTH: 3451

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 18

```

tcgagacggg aagccactct ctacgagaag acagaagccc ctcacagagg cctctgtcta    60

```

-continued

cgccactaa agctcggctt attcatatgt atttatattc tttcaataga tcaactcagcg	120
ctatthtaag ttcacctctt gtaagtccac ctgggcgctc tttctttcct tcggtaaagc	180
tgctggccag accaaacatt aaactcaagc atctcccaag cgatgcatca tcttgggcca	240
gcatccctga atcgcgcgct ggacctccaa gtcttaaaaa attcttcgct gaaggttttc	300
ccatcaatcg atgaggctaa tagcttcttt gcaatatcta tcatttccat gctcacctta	360
aagcacctca tttttcatgt aaaaattgta ttgatccgtg ccagactcaa tcctccacc	420
agaacaaca atcccatcct ctccaatgat aacaacaata ttagtccctgg cattgtaatg	480
tacttttgag ttacttcgag agtggttaagt ccctttttct acggttgcaag gatcagcaag	540
gtgctcaaga attttatccc taaactctgc aagcgttcca ttgttggcgc ttttttcacc	600
cagcccaaaa tcatatthgt ggctatcaaa ttttttctgt agttgcctcc gttgtaagat	660
accactatca agaggactac tgagcattac ataaacaggt ttgactccag aatccgccgg	720
gaaaatcacg atcagatcgt ttaggccag tagcattccc ggataggact ccgggccggg	780
cttcaacggt gtgagggccg ctccctcata taccggcacc ggcttcggta tgaccggagt	840
ggtactcga gggttctggt ttcctggagg actcgcggc gtccaagtca ggatcagtgg	900
cggcgcttct gcgaccgtag agggaaccgt aacctcgtac agtccctgtg cggcgttata	960
ggccccatcc ggaccggaac gctttcggaa cgctcacacc atcggctga ccaccgaaag	1020
gtcgtcgtgt tgctcgcgc ctctgtggtc aggcgcctcg gcagatcgac ggtaccgctg	1080
gcttttgcaa ccgcttcag gtttacgctt gggggaagcc ccaatttagc ggcattcatg	1140
cccagggcgt aacgaacgct atcggggcgt ttggtcctgcc attgctcggc agtccgggag	1200
agtaggtcag actggcaagc cacggccatc accgaggtgc tgaagccagg accgccagga	1260
cggcaatcgc atcggagatc gcttgagcaa gggatgcggc gcctgtgcga cctggatcag	1320
accccgctgc ggcgggtggc caccgcctgc cattggctgg catggcataa gtattggcag	1380
ccctgatcgc cgettgacga gcgatttcct tgcgccttgc cgtttcggcg ttcagcttgt	1440
ccagccgtgc ttgcaggctg gcgatttcat ccactagta ggacatcggc gttgtaggtt	1500
gccttttggt tctccagtgc attgggtgcc ttggcaatca aggcattggt tgcagtctgc	1560
aattcttctt attgcgatcg cctgcgtaag gagttgagta gcgcgttcaa gccactgctc	1620
tggcgttga ttggtcagtt gaggcaaac attcccagcc tggtaagct cggactgcac	1680
ttttttctcg acatttgctt tcttggcctt gtagtccgcc tccacctcag cagcggctcg	1740
ctgggcttct gcttccaatg accgggcttt attctccagc tcttgagacg tttgthtcaa	1800
gatagcgatt tgcgccttat agatcggc gctgtacgct ttggccagct cactcatatg	1860
gcgatccag aactctccat agaattttcg gctggccagc aactgactct ggtacatcga	1920
ctctgacttc tgaggaaagt ctgaagccgt ataaagattg gccggcgat cctcaatgac	1980
ctttagcgat ttgtctttg catccatgag tgcatacaag atactctttt catcgcggat	2040
gtcattggca ctgaccgctt tacctggcaa ccccgcttca ctcttgagtt catcaacctc	2100
cttcagggtt tcatttttca ggtttttctt gagttctgaa tgggacttat caagcgtact	2160
tcttagcttc ctgtactcct gcattocagt accgacatac ggacttggtc ctggtgggac	2220
aaatggtgga gtaccgtagc ttgatcagc aggaatatac tggattatgt cacgccacc	2280
accctgcaca tgtgtaataa ccatacgaacc aggttcgtaa tcattgacag ccatagatcg	2340

-continued

```

ccccacatt aatttgaaag tgtaatgtat tgagcgactc ccacctagag aacctctctc 2400
cagtcaataa gcccacatgc atcggcaata cactgcaatc aacttcaata tcccggtgtt 2460
agatgatcca gaaggtgcgc tctctcgcct cttataatcg cgctcgctgc aaacggtcac 2520
ttccttaacg cacacctcat ctaccccggc cagtcacgga agccgcatac cttcggttca 2580
ttaacgaact cccactttca aaattcatcc atgccgcccc ttcgagagct tccggacaaa 2640
gccacgctga ttgagagccc agcgtttttg attgcaagcc gctgcagctg gtcaggccgt 2700
ttccgcaacg cttgaagtcc tggccgatat accggcaggg ccagccatcg ttcgacgaat 2760
aaagccacct cagccatgat gccctttcca tccccagcgg aaccccgaca tggacgccaa 2820
agccctgctc ctccggagcc tctgectggc cgccccattc gccgacgcgg cgacgctcga 2880
caatgctctc tccgctgcc tcgccgcccg gctcgggtgca ccgcacacgg cggagggcca 2940
gttgcacctg cactcaccct ttgagggccc gcgctccacc gggaatgag gctgtacctc 3000
ggcgctgggt cgatatcgcc tgctggccag gggcgccagc gccgacagcc tegtgttca 3060
agagggctgc tcgatagtcg ccaggacacg ccgcgcacgc tgaccctggc ggcggacgcc 3120
ggcttggcga gcggccgcga actggctcgc accctgggtt gtcaggcgcc tgactgacag 3180
gccgggctgc caccaccagc ccgagatgga cgccctgcat gtatcctccg atcggcaagc 3240
ctcccgcttc cacattcacc actctgcaat ccagttcata aatcccataa aagccctctt 3300
ccgctccccg ccagcctccc cgcctccccc accctagacg ccccgccgct ctccgcccgc 3360
tcgcccgaca agaaaaacca accgctcgat cagcctcatc cttaccatcac cacaggagcc 3420
atcgcgatgc acctgatacc ccattggatc c 3451

```

```

<210> SEQ ID NO 19
<211> LENGTH: 744
<212> TYPE: DNA
<213> ORGANISM: Pseudomonas aeruginosa

```

```

<400> SEQUENCE: 19

```

```

gggttcagca agcgttcagg ggcggttcag taccctgtcc gtactctgca agcgtgaaac 60
gacacgactc tcgcagaacg gagaaacacc atgaaagcac tcaagactct ctctcatgcc 120
accgccctgc tgggttccgc cgccggcgctc caggccgcgg acaacttctg cgccctgacc 180
tggggcgaga ccagcaacaa catccagaaa tccaagtgcg tgaaccgcaa cctgaacagc 240
ccgaacctcg acaaggtgat cgacaacacc ggcacctggg gcctccgccc cggccagcag 300
ttcgagcagg gccgctacta cgcgacctac gagaacatct ccgacaccag cagcggcaac 360
aagctgcgcc agcagaacct gctcggcagc tacgacgctt tcctgccgat cggcgacaac 420
aacaccaagc tgctcggcgg tgccaacctc ggcctggtca agctggaaca ggaocggcaag 480
ggcttcaagc gcgacagcga tctcggctac gctgcggggc tgacggccgg tctcctgcag 540
gagctgagca agaatgcctc gatcgaaggc ggctatcgtt acctgcgcac caacgccagc 600
accgagatga ccccgcatgg cggcaacaag ctgggctccc tggacctgca cagcagctcg 660
caattctacc tgggcgccaa ctacaagttc taaatgaccg cgacgcgccc gcgagggcat 720
gcttcgatgg ccgggcccga aggt 744

```

```

<210> SEQ ID NO 20
<211> LENGTH: 2760
<212> TYPE: DNA

```

-continued

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 20

```
ctgcagctgg tcaggccgtt tccgcaacgc ttgaagtccct ggccgatata ccggcagggc      60
cagccatcgt tcgacgaata aagccacctc agccatgatg ccctttccat ccccagcggga    120
accccgacat ggacgcaaaa gccctgctcc tcggcagcct ctgcctggcc gccccattcg    180
ccgacgcggc gacgctcgac aatgctctct ccgctgcct cgcgcccggg ctcggtgcac    240
cgcacacggc ggagggccag ttgcaactgc cactcaccct tgaggcccgg cgtccaccg     300
gcgaatgcgg ctgtacctcg gcgctggcgc gatatcgct gctggccagg ggcgccagcg    360
ccgacagcct cgtgcttcaa gagggctgct cgatagtcgc caggacacgc cgcgcacgct    420
gaccctggcg gcggacgccg gcttggcgag cggccgcgaa ctggtcgtca cctggggttg    480
tcaggcgcct gactgacagg ccgggctgcc accaccaggc cgagatggac gccctgcatg    540
tatactccga tcggcaagcc tcccgttcgc acattcacca ctctgcaatc cagttcataa    600
atcccataaa agccctcttc cgtccccgc cagcctcccc gcatacccga ccctagacgc    660
cccgcgcctc tccgcgggct cgcgcgacaa gaaaaaccaa ccgctcgatc agcctcatcc    720
ttcaccatc acaggagcca tcgcatgca cctgataccc cattggatcc cctggctcgc    780
cagcctcggc ctgctcggc gcggctcgtc cgcgtccgcc gccgaggaag ccttcgacct    840
ctggaacgaa tcgcccctaa cctgcgtgct cgacctcaag gacggcgtgc gttccagccg    900
catgagcgtc gacccggcca tcgcccacac caacggccag ggcgtgctgc actactccat    960
ggtcctggag ggcggcaacg acgctcctca gctggccatc gacaacgccc tcagcatcac   1020
cagcgcggc ctgacctcc gcctcgaagg cggcgtcgag ccgaacaagc cggcgcgcta   1080
cagctacacg cgcaggcgc gcggcagttg gtcgctgaac tggctggtac cgatcggcca   1140
cgagaagccc tcgaacatca aggtgttcat ccacgaactg aacccggca accagctcag   1200
ccacatgctc ccgatctaca ccacgagat gggcgacgag ttgctggcga agctggcgcg   1260
cgatgccacc ttctctgca gggcgacga gagcaacgag atcgagccga cgtcgcctat   1320
cagccatgcc ggggtcagcg tggctatggc ccagaccag ccgcccggg aaaagcgtg    1380
gagcgaatgg gccagcggca aggtgtttgt cctgctcgac ccgctggacg gggcttacia   1440
ctacctgcc cagcaacgct gcaacctcga cgatacctgg gaaggcaaga tctaccgggt   1500
gctgcggcgc aacccggcga agcatgacct ggacatcaaa cccacggcga tcagtcatcg   1560
cctgcacttt cccgagggcg gcagcctggc cgcgctgacc gcgcaccagg cttgccacct   1620
gccgctggag actttcacc gtcatcgcca gccgcggcgc tgggaacaac tggagcagtg   1680
cggctatccg gtgcagcggc tggctgcct ctacctggcg gcgcggtgt cgtggaacca   1740
ggtcgaccag gtgatccga acgcccctgc cagccccggc agcggcggcg acctgggcga   1800
agcgtaccg gagcagccgg agcagggccg tctggccctg acctggccg ccgcccagag   1860
cgagcgttc gtccggcagg gcaccggcaa cgacgaggcc ggcgcggcca acgcccagct   1920
ggtgagcctg acctgcccgg tcgcccggcg tgaatgcgcg ggcggcggcg acagcggcga   1980
cgccctgctg gagcgcgaact atcccactgg cgcggagtcc ctcggcgacg ggcggcagct   2040
cagcttcagc acccgcggca cgcagaactg gacggtggag cggctgctcc aggcgcaccg   2100
ccaactggag gagcgcggct atgtgttcgt cggctaccac ggcacctcc tcgaagcggc   2160
```

-continued

```

gcaaagcacc gtcttcggcg ggggtcgcgc gcgcagccag gacctcgacg cgatctggcg 2220
cggtttctat atcgccggcg atccggcgct ggcctacggc tacgcccagg accaggaacc 2280
cgacgcacgc ggccggatcc gcaacgggtgc cctgctgcgg gtctatgtgc cgcgctcgag 2340
cctgcccggc ttctaccgca ccagcctgac cctggccggc ccggaggcgg cgggcgaggt 2400
cgaacggctg atcgccacc cgctgcccct gcgcctggac gccatcaccg gccccgagga 2460
ggaaggcggg cgcctggaga ccattctcgg ctggcccgct gccgagcgca ccgtggtgat 2520
tcccctcgcg atccccaccg acccgcgcaa cgtcggcggc gacctcgacc cgtccagcat 2580
ccccgacaag gaacaggcga tcagcgcctt gccggactac gccagccagc ccggcaaacc 2640
gccgcgcgag gacctgaagt aactgcccg accggccggc tcccttcgca ggagccggcc 2700
ttctcggggc ctggccatac atcaggtttt cctgatgcca gcccaatcga atatgaattc 2760

```

```

<210> SEQ ID NO 21
<211> LENGTH: 172
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

```

```
<400> SEQUENCE: 21
```

```

ttgatgaaat gcatcgatta ataaattttc atgtacgatt aaaacgtttt tacccttacc 60
ttttcgtact acctctgcct gaagttgacc acctttaaag tgattcgttg aaatccatta 120
tgctcattat taatacagtc tataaaaaca aatggaatgt gatgatcgat ga 172

```

```

<210> SEQ ID NO 22
<211> LENGTH: 155
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

```

```
<400> SEQUENCE: 22
```

```

gttcattga ctctgtatca cctgttgtaa cgaacatcca tatgtcctga aactccaacc 60
acaggtttga ccacttccaa tttcagacca ccaagtttga cacgtgaaga ttcattcttct 120
aatatttcgg aattaatatc atattattta aatag 155

```

```

<210> SEQ ID NO 23
<211> LENGTH: 145
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

```

```
<400> SEQUENCE: 23
```

```

acatagaaa actcaaaaga tttacttttt tcaaatggaa aataagggta cacacgatat 60
ttcccgatc cttcagttac cggtaacaac tctcttttat taacctgcac ataactctgac 120
tccgcttcac tcatcaaact actaa 145

```

```

<210> SEQ ID NO 24
<211> LENGTH: 266
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

```

```
<400> SEQUENCE: 24
```

```

tttcaactga attacatttc gctcattacg tacagtgaca atcgcgctcag atagtttctt 60
ctggttagct tgactcttaa caatctgttc taaattttgt ttaattcttt gattcgtact 120
agaaatttta cttctaattc cttgtaattc ataacttgca ttatcatata aatcataagt 180

```

-continued

```
atcacatddd tgatgaatac tttgatataa atctgacaat acaggcagtt gctccattct 240
atcgttaaga atagggtaat taatag 266
```

```
<210> SEQ ID NO 25
<211> LENGTH: 845
<212> TYPE: DNA
<213> ORGANISM: Haemophilus influenzae
```

```
<400> SEQUENCE: 25
```

```
tgttaaattd ctttaacagc gattttgtda tttaaattaa acctattatt ttgtcgcttc 60
tttcaactgc tctactgctt gagttgcttt ttctgaaacc gcctctttca tttcaactgc 120
ttttctgat gctgcttctt tcaattcgcc tactttttct gacgctgctt ctggtgctga 180
tttaattact tctttcgcac cttccacttt ctctgctact ttatttttca cgtctgtaga 240
aagctgctgt gctttttcct ttacttcagt cattgtatta gctgcagcat cttttgtttc 300
tgatgcgact gatgctacag ttgcttcgt atcctcaact tttgtttttg cttcttgctt 360
atcaaaaaca cctgtcacga ctaaagctga acctaaaacc aatgctaag ttaatttttt 420
cattattdtc tccatagaat aatttgattg ttacaagcc ctattacttt gatgcagttt 480
agtttacggg aattttcata aaaagaaaa cagtaatagt aaaactttac ctttctttaa 540
aaagattact ttataaaaa acatctaaga tattgatttt taatagatta taaaaacca 600
ataaaaattd tattttttgt aaaaaaaag aatagtttat tttaaataaa ttacaggaga 660
tgcttgatgc atcaatattt ctgatttatt accatcccat aataattgag caatagttgc 720
aggataaaat gatattggat ttctgtttcc atacagttca gcaacaattd ctcccactaa 780
gggcaaatgg gaaacaatta atacagattd aacgccctcg tcttttagca cttctaaata 840
atcaa 845
```

```
<210> SEQ ID NO 26
<211> LENGTH: 1598
<212> TYPE: DNA
<213> ORGANISM: Haemophilus influenzae
```

```
<400> SEQUENCE: 26
```

```
gaatagagtd gcaactcaata gattcgggct ttataattgc ccagattdtd atttataaca 60
aagggttcca aatgaaaaa tttaatcaat ctctattagc aactgcaatg ttgttgctg 120
cagggtgtdc aaatgcggca gcggttcaat tggcggagtd ttctacttca ggtcttgtdc 180
gtgcctatgc ggggtgaagc gcgattgcag ataatgcttc tctctgtdgc actaaccag 240
ctttgatgag tttattttaa acggcacagtd tttccacagc tggcgtttat attgattcta 300
gaatataat gaatggtgat gtaacttctt atgctcagat aataacaaat cagattggaa 360
tgaaagcaat aaaggacgdc tcagcttcac agcgtaatgt tgttcccggd gcttttgtdc 420
caaatcttda tttcgttgdc ccagtgatg ataaattgdc gctgggtgct ggaatgaatg 480
tcaattdcgc tctaaaaagt gaatatgacg atagttatga tgctggtgta tttggtggaa 540
aaactgactt gagtgcatac aacttaaat taagtgtgdc ttatcgagta acagaaggtt 600
tgagcctagc ttttagggta aatgcggttd atgctaaagc ccaagttgaa cggatgctg 660
gtcttattgc ggatagtdt aaggataacc aaataacaag cgcactctca acacagcaag 720
aaccattcag agatcttaag aagtattgdc cctctaagga caaatctgtd gtgtcattac 780
```

-continued

```

aagatagagc cgcttggggc tttggctgga atgcagggtg aatgtatcaa tttaatgaag 840
ctaacagaat tggtttagcc tatcattcta aagtggacat tgattttgct gaccgcactg 900
ctactagttt agaagcaaat gtcacaaaag aaggtaaaaa aggtaattta acctttacat 960
tgccagatta cttagaactt tctggtttcc atcaattaac tgacaaaact gcagtgcaat 1020
atagttataa atatacccat tggagtcggt taacaaaatt acatgccagc ttcgaagatg 1080
gtaaaaaagc ttttgataaa gaattacaat acagtaataa ctctcgtggt gcattagggg 1140
caagttataa tctttatgaa aaattgacct tacgtgctgg tattgcttac gatcaagcgg 1200
catctcgtca tcaccgtagt gctgcaattc cagataccga tcgcacttgg tatagtttag 1260
gtgcaaccta taaattcacg ccgaatttat ctggtgatct tggctatgct tacttaaaag 1320
gcaaaaaagt tcactttaaa gaagtaaaaa caataggtga caaacgtaca ttgacattga 1380
atacaactgc aaattatact tctcaagcac acgcaaatct ttacggtttg aatttaaatt 1440
atagtttcta atccgttaaa aaatttagca taataaagca caattccaca ctaagtgtgc 1500
ttttctttta taaaacaagc cgaaaaatga ccgcacttta ttacacttat taccctcgc 1560
cagtcggacg gcttttgatt ttatctgacg gcgaaaca 1598

```

<210> SEQ ID NO 27

<211> LENGTH: 9100

<212> TYPE: DNA

<213> ORGANISM: Haemophilus influenzae

<400> SEQUENCE: 27

```

gtcaaaaatt gcgtgcattc tagcgaaaaa atgggctttt gggaaactgtg ggatttattt 60
aaaatcttag aaaatcttac cgcactttta agctataaag tgcggtgaaa tttagtggcg 120
ttataaatgg agaattactc tgggtgaatc cattcgactg tccagcttcc agtaccttct 180
ggaaactaat tttttgtgag ataaggcaaa atttctttca tttgggtttc taatgtccaa 240
ggtggattaa ttaccacat accgctcgca gtcattcctc gttgatcgtc atctgggcga 300
acggcgagtt caatttttag aatttttcta attcccgttg cttctaaacc cttaaaaata 360
cgttttagtt gttggcgtaa tacaaacgga taccaaatcg cataagtgcc agtggcaaaa 420
cgtttatagc cctcttcaat ggctttaaca acgagatcat aatcatcttt taattcataa 480
ggcggatcga tgagtactaa gcctcggcgt tcttttgcg gaagcgttgc tttgacttgt 540
tgaaagccat tgtcacattt tacggtgaca tttttgtcgt cgctaaaatt attgcaaga 600
attggataat cgctaggatg aagctcggtc aatagtgcgc gatcttgtga gcgcaacaat 660
tccgcggcaa ttaatggaga acccgcgtaa taacgtagtt ctttgccacc ataattgagt 720
tttttgatca tttttacata acgagcaata tcttcgggta aatctgtttg atcccacagg 780
cgtccaatac cttctttata ttccccgtt ttttctgatt catttgagga taaacgataa 840
cgccccacac cagagtgcgt atccaaaata aaaaagcett tttctttgag ttttaagattt 900
tccaaaaatga gcattaaaac aatatgtttc aagacatcgg catgattgcc agcgtgaaat 960
gagtgatgat aactcagcat aatatattcc ttatatattc cttatttgtt taataacgaa 1020
ggcggacca ttgactcggc cgattacaca ctaaagtgcg gtcattttta gaagagttct 1080
tgtggttgcg tcgctggcgt attgcottca ttatttaagc gttgctgtaa ctcagtagga 1140
acataataac cacgctcttg catttccgaa agataggtac gtgtcggttc tgttcccgca 1200

```

-continued

ataaaatatt	ctttgcgcc	accgtttga	gaaagcaaac	ctgtcaaagt	atcaatgttt	1260
ttttccaaa	tttttgccg	tagcgacaat	ttacgttctg	gcttatcact	caaagccgtt	1320
ttcatataag	tgatccaagc	aggcattgct	gtttttgctc	ctgcttctcc	acgcccaagt	1380
actggtttgt	tatcatcaaa	cccgaacata	gttggtggtta	ctaagtttgc	accaaattccc	1440
gcataccaag	coacttttga	actgttggtta	gtacctgttt	taccgcctat	atcgctacgt	1500
ttaatgcttt	gtgcaatacg	ccagctggtg	cctttccagt	ctaaaccttg	ttcgccataa	1560
attgccgtat	ttaaaggcact	acgaatgaga	aaagcaagtt	cgccactaat	gacacgtggc	1620
gcataattcta	ttttcgacga	agcatttttt	gcagcagcca	ttaaatcaat	cgcatcttct	1680
ttaagtgcgg	tcatatttga	ttgtaattct	ggcagttcag	gcacagtttc	aggttgttga	1740
tctaattctt	cgccattggt	gctgtcatct	gttggtttta	aggcattctc	gcctaaagga	1800
atattggcaa	agccgttgat	tttgtctttg	gtttcgccat	aaattacagg	tatatcatta	1860
cattcaatgc	aagcaatttt	agggtttgca	ataaataagt	ctttaccctg	gttatcttga	1920
atTTTTTcaa	tgatataag	ttcaatgagg	aagccaccat	tatcaaacac	cgcataagct	1980
cgcgccattt	ctaagtgtgt	gaaagaggct	gcgccaagt	ctaaggcttc	actggcaaaa	2040
tattgatcac	gtttaaaacc	aaaacgttgt	aaaaattctg	ctgtgaaatc	aatacctgcc	2100
gtttgatag	cacgaatagc	aattatattt	ttggattgac	ctaactctac	gcgtaaaccg	2160
atcgggccat	cataacgatc	aggcgagttt	ttcggttgcc	acattttttg	tcccggtttt	2220
tgaatagaaa	tcgggctgtc	ttgtaatacg	cttgaaagt	ttaagccttt	ttctaattgt	2280
gccgcgtaaa	taaatggttt	gatagaagaa	cccacttgaa	ctaaagactg	tgtggctcga	2340
ttgaatttac	ttgttcata	gctaaagcca	ccgaccactg	cttcaatcgc	accattatct	2400
gaattaagag	aaactaatgc	tgaatttgc	gcgggaattt	gtcctaattg	ccattcccca	2460
ttagcacgct	gatgaatcca	aatttgcctg	ccgactttca	caggattgct	tctgcctgtc	2520
caacgcattg	cattggttga	taaggtcatt	ttttcccag	aagcgagcaa	tatatcagca	2580
ccgcctttta	caattccaat	cactgccgca	ggaataaatg	gctctgaatc	aggtagtgtg	2640
cgtagaaaa	cgacaatg	atcattgtcc	caagcggctt	catttttttg	ccataatggc	2700
gcgccaccgc	gataaccgtg	acgcatatcg	taatcaatca	agttattacg	cacagctttt	2760
tgggcttcag	cttggtcttt	tgaaagtaca	gtggtaaata	ctttataacc	actggtgtaa	2820
gcattttctt	cgccaaaacg	acgcaccatt	tottgacgca	ccatttcagt	gacataatcg	2880
gctcgaaatt	caaattttgc	gccgtgatag	ctcgcccaa	tcggctcttt	caatgcagca	2940
tcatattctt	ctttgctgat	gtatttttca	tctaacatac	ggcttagcac	cacattgccc	3000
cgttcttctg	aacgttttaa	agaataaagc	gggttcattg	ttgaaggtgc	tttagttaa	3060
ccagcaataa	tcgccatttc	cgataaggtc	aattcattca	atgatttacc	gaaatagggt	3120
tgtgctgccg	ctgcaacacc	ataagaacga	tagcctaaaa	agattttggt	taaataaagc	3180
tctaataattt	cttggttgtt	gagagtattt	togatttcta	ccgcaagcac	ggcttcacga	3240
gctttacgaa	taatggtttt	ttctgaggtt	aagaaaaagt	tacgcgctaa	ttgttgagta	3300
atcgtaactg	cgcttctgta	tgcaccgcca	ttactcactg	cgcaaaaaa	tgcacgggca	3360
atgccgatag	ggtctaattc	gtgatgatcg	taaaaacgac	tgtcttcogt	cgtaaaaaat	3420
gcgtcaatta	agcgttggg	cacatcggt	aatttcactg	gaatacggcg	ttgctcacc	3480

-continued

acttcgcca	ttaatttacc	gtcagccgta	taaatctgca	ttggttgctg	taattcaacg	3540
gtttttaatg	tttctactga	gggcaattca	gattttaagt	ggaatacaa	cattccgcct	3600
gctactaaac	ctaaaataca	taaagttaat	agggtgttta	atattaatth	tgcgatccgc	3660
atcgtaaaat	tctcgtctcg	ttaatgaata	ttcttgtcaa	gagacctatg	atttgctgt	3720
taagtataaa	agattcagcc	tttaagaat	aggaaagaat	atgcaattct	ccctgaaaa	3780
ttaccgcact	ttacaaatcg	gcattcatcg	taagcagagt	tattttgatt	ttgtgtggtt	3840
tgatgatctc	gaacagccac	aaagtatca	aatctttggt	aatgatcgtt	atthtaaaaa	3900
tcgtttttta	caacagctaa	aaacacaata	tcaagggaaa	acctttcctt	tgcagtttgt	3960
agcaagcatt	cccgcccact	taacttggtc	gaaagtatta	atgttgccac	aagtgttaa	4020
tgcgcaagaa	tgatcatcaac	aatgtaaatt	tggtgattgaa	aaagagctgc	ctatthtttt	4080
agaagaattg	tggtttgatt	atcgttctac	cccgttaaag	caaggthttc	gatttagagg	4140
tactgcaatt	cgtaaaagta	gcgctcaaac	ttatthgcaa	gattthtcagc	catttaatat	4200
taatatattg	gatgttgcgt	caaagtctgt	tttgctgca	thtcaatath	tgthgatga	4260
acaagtgcgg	tcagaaaata	ccttatthtt	atthcaagaa	gatgactatt	gcttgccgat	4320
ttgtgaaaga	tctcagcaat	cacaaattht	acaatctcac	gaaaatthga	ccgcacttht	4380
tgaacaatth	accgaacgth	ttgaaggaca	acttgaacaa	gthtttgtht	atcaaatthc	4440
ctcaagtcac	acaccattac	ccgaaaactg	gcagcgagta	gaaacagaac	thcctthtat	4500
tgcgctgggc	aacgcgctat	ggcaaaaaga	thtathcatca	caaaaagthg	gtgthtaaat	4560
gtcgtatgat	thattgcctt	ggcgtactth	tcaacatca	aagcgtthac	gtcgtthtag	4620
thtttatath	gctthattth	tcttgccttg	tattaattht	atgthgctt	thagcaatth	4680
gattgaacaa	cagaaacaaa	atthtgcagc	acagcaaaag	tcgthtgaac	aacttaatca	4740
acagcttcat	aaaactacca	tgcaaatthg	tcagthtacgc	atthcggthg	aagthgthg	4800
agthttgaca	tctatthcca	acgagcaagt	aaaaagagth	thacaacagc	taagthgaat	4860
acctthtcaa	caaggagaac	tgaataaatt	taaaacagat	gccaataact	taagctthg	4920
aggtaacgcg	caagatcaaa	cagaatthg	actgattcat	caatthttaa	agaacatth	4980
thccaatgth	aaatthagth	aggthcaacc	tgaacaagat	acatthgtht	thcactthg	5040
tgthgaacaa	gggcccggaa	aatgaaagct	thttthtaacg	atcctthtac	thcctthg	5100
aaatggctaa	gtcagccctth	thatgtgcac	ggtthaacct	thttatthgct	atthagthgcg	5160
gtgaththtc	gccccgthtth	agattatata	gaggggagth	cacgththcca	tgaathgaa	5220
aatgagthtag	cggtgaaacg	thcagaatthg	thgcatcaac	agaaaattht	aacctcttht	5280
caacagcagth	cgaaaagthc	aaaactthct	ccagaactgg	ctgcacaaat	tatthcctthg	5340
aataaacaaa	thcaacgthtth	agctgcgcgth	aacgththtth	ctcagcatht	acgthgggaa	5400
atggggcaaa	agcctatthtth	gcattthacag	ctthacagthc	atthtthgaaa	aacgaagaca	5460
thtttatccg	cactthtthg	taattogthca	cagctthctg	taagthcggth	gcaatthtag	5520
aaacccgaag	acggcccatt	gcaaaccgag	atcaththtthc	agctagataa	ggaacaaaa	5580
tgaacatthg	gthttthcctg	atthattatth	thttthtagaa	thgagthtthg	ggacaagathc	5640
ctthcagataa	aacacagcgt	aaccgthctc	agthtthgataa	cgcaacaaaca	gtaatggagc	5700
aaacagaaat	aatthctctca	gatgtgccta	ataatctatg	cgagcggat	gaaaatcgcc	5760

-continued

aagcggctga aattcctttg aacgctttaa aattggtggg ggtagtgatt tctaaagata	5820
aagcctttgc ctgtttgcaa gatcaaggtt tgcaagttaa cagcgtttaa gagggcgttg	5880
atgtggctca agaggctat attgtagaaa aaatcaacca aaacaatggt caatttatgc	5940
gtaagctagg agagcaatgt gatagtagtg aatggaaaaa attaagtttt taaaggaaga	6000
ttatgaagaa atatTTTTTt aagtgcggtt atTTTTtagt atgTTTTgt ttgccattaa	6060
tcgTTTTtgc taatcctaaa acagataacg aacgTTTTtTt tattcgttta tcgcaagcac	6120
ctttagctca aacactggag caattagctt ttcaacaaga tgtgaattta gtgattggag	6180
atataattgga aaacaagatc tctttgaaat taaacaatat tgatatgcca cgtttgctac	6240
aaataatcgc aaaaagtaag catcttactt tgaataaaga tgatgggatt tattatttaa	6300
acggcagtca atctggcaaa ggtcaggttg caggaaatct tacgacaaat gaaccgcact	6360
tagtgagtca cacggtaaaa ctccattttg ctaaagcttc tgaattaatg aaatccttaa	6420
caacaggaag tggctctttg ctttctcccg ctgggagcat tacctttgat gatcgagta	6480
atttgctggt tattcaggat gaacctcgtt ctgtgcaaaa tatcaaaaaa ctgattgctg	6540
aaatggataa gcctattgaa cagatcgcta ttgaagcgcg aattgtgaca attacggatg	6600
agagtttgaa agaacttggc gttcgggtgg ggatttttaa tccaactgaa aatgcaagac	6660
gagttgcggg cagccttaca ggcaatagct ttgaaaaat tgcggataat cttaatgtaa	6720
atTTTtgcgac aacgacgaca cctgctggct ctatagcatt acaagtcgcc aaaattaatg	6780
ggcgattgct tgatttagaa ttgagtcgct tggagcgtga aaataatgta gaaattattg	6840
caagccctcg ctactcact accaataaga aaagtgcgag cattaacag gggacagaaa	6900
ttccttacet cgtgagtaat actcgtaacg atacgcaatc tgtggaattt cgtgaggcgg	6960
tgcttggttt ggaagtgcg ccacatattt ctaaagataa caatatctta cttgatttat	7020
tggtaaagca aaattcccct ggttctcgtg tcgcttatgg acaaaatgag gtggtttcta	7080
ttgataaaca agaaattaat actcaggttt ttgccaaaaga tggggaacc attgtgcttg	7140
gcggcgattt tcacgataca atcacgaaaa gcgaagataa agtgccattg cttggcgata	7200
taccgcttat taaacgatta tttagcaaaag aaagtgaacg acatcaaaaa cgtgagctag	7260
tgattttcgt cacgccacat atTTTttaaag caggagaaaa cgttagaggc gttgaaacaa	7320
aaaagtgcg gtaaaaaata actTTTTttaa tgatgaattt tTTtatttt cgtgtatcc	7380
actgtcgtgg caatcttcat atcgcaaaaa atgggttatg ttcaggttg caaaaacaaa	7440
ttaaatcttt tccttattgc ggtcattgtg gttcgggaatt gcaatattat gcgcagcatt	7500
gtgggaattg tcttaacaaa gaaccaagtt gggataagat ggtcattatt gggcattata	7560
ttgaacctct ttcgatattg attcagcgtt ttaaatttca aaatcaattt tggattgacc	7620
gcactttagc tcggctttta tatcttgcg tacgtgatgc taaacgaacg catcaactta	7680
aattgcaga ggcaatcatt ccagtcctt tataatcattt tcgtcagtg cgacggggtt	7740
ataatcagc agatttatta tctcagcaat taagtcgttg gctggatatt cctaatttga	7800
acaatatcgt aaagcgtgtg aaacacacct atactcaacg tggtttgagt gcaaaagatc	7860
gtcgtcagaa tttaaaaaat gccttttctc ttgctgtttc gaaaaatgaa tttccttate	7920
gtcgtgttg gttggtggat gatgtgatta ctactggttc tactactcaat gaaatctcaa	7980
aattgttgcg aaaattaggt gtggaggaga ttcaagtgtg gggctggca cgagcttaat	8040

-continued

```

ataaagcact ggaaaaaaaa gcgcgataag cgtattatcc cgcatacttt ctctcaagta 8100
tttaggacat aattatggaa caagcaaccc agcaaatcgc ttttctgat gccgcacaag 8160
cgcattttcg aaaactttta gacacccaag aagaaggaac gcatattcgt attttcgcgg 8220
ttaatcctgg tacgcctaata gcggaatgtg gcgtatctta ttgccccccg aatgccgtgg 8280
aagaaagcga tattgaaatg aaatataata ctttttctgc atttattgat gaagtgagtt 8340
tgcccttctt agaagaagca gaaattgatt atgttaccga agagcttggg gcgcaactga 8400
ccttaaaagc accgaatgcc aaaatgcgta aggtggctga tgatgcgcca ttgattgaac 8460
gtgttgaaata tgtaattcaa actcaaatta acccacagct tgcaaatcac ggtggacgta 8520
taaccttaat tgaattact gaagatggtt acgcagtttt acaatttggg ggtggctgta 8580
acggttggtc aatggtggat gttacgtaa aagatggggg agaaaaacaa cttgttagct 8640
tattcccgaa tgaattaaaa ggtgcaaaag atataactga gcatcaacgt gccgaacatt 8700
cttattatta gtgagttata aaagaagatt tataatgacc gcacttttga aagtgcggtt 8760
atttttatgg agaaaaaatg aaaatacttc aacaagatga ttttggttat tggttgctta 8820
cacaaggttc taatctgtat ttagtgaata atgaattgcc ttttggtatc gctaaagata 8880
ttgatttga aggattgcag gcaatgcaaa ttggggaatg gaaaaattat ccgttgggc 8940
ttgtggctga gcaagaaagt gatgaacgag aatatgtgag tttgagtaac ttgctttcac 9000
tgccagagga tgaattccat atattaagcc gaggtgtgga aattaatcat tttctgaaaa 9060
cccataaatt ctgtggaaag tgcggtcata aaacacaaca 9100

```

```

<210> SEQ ID NO 28
<211> LENGTH: 525
<212> TYPE: DNA
<213> ORGANISM: Moraxella catarrhalis

```

<400> SEQUENCE: 28

```

aaaaatcgac tgccgtcatt ttcaaccacc acatagctca tattcgcaag ccaatgtatt 60
gaccgttggg aataataaca gccccaaaac aatgaaacat atggtgatga gccaaacata 120
ctttcctgca gattttggaa tcatatgcc atcagcacca gtatggtttg accagtattt 180
aacgccatag acatgtgtaa aaaaattaaa taacggtgca agcatgagac caacggcacc 240
tgatgtacct tgtacgatga cctcacctgc tgtggcaacc ataccaagtc cattgcctgt 300
gatatttttg cgaaaagaca aacttaccac acagaccaag ccgatgattg agatgacaaa 360
ataaaaccaa tccaaatgcg tgtgagctgt tgtggtccaa aatccagtaa atagtgcaat 420
aaatccgcaa acaaaacaaa gtagcaccca gcttgttgc caatcttttt taccaaagcc 480
tgtgatgtta tctaaaatat caattttcat cagattttcc ctaat 525

```

```

<210> SEQ ID NO 29
<211> LENGTH: 466
<212> TYPE: DNA
<213> ORGANISM: Moraxella catarrhalis

```

<400> SEQUENCE: 29

```

taatgataac cagtcaagca agctcaaatc agggtcagcc tgttttgagc tttttatttt 60
ttgatcatca tgcttaagat tcaactctgc atttttttac aacctgcacc acaagtcatc 120
atcgcatthy caaaaatggt acaaaacaag cgtcagcgac ttaaacaaaa aaaggctcaa 180

```

-continued

tctgctgtg tgcgttcaact ttacaaaac accatgcacc gctttgacat tgttggtgaa	240
tttcatgacc atgcacaccc ttattatatt aactcaaata aaatacgeta ctttgtcagc	300
tttagccatt cagataatca agtcgctctc atcatcagct taacaccttg tgccattgac	360
atagaagtta acgatattaa atacagtgtg gttgaacgat actttcatcc caatgaaatt	420
tatctactta ctcaatttag ctctactgat aggcaacagc ttatta	466

<210> SEQ ID NO 30
 <211> LENGTH: 631
 <212> TYPE: DNA
 <213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 30

gatctttgat tttcatgtag tattactctc tcttgtcact tctttctatt ttaccataaa	60
gtccagcctt tgaagaactt ttactagaag acaaggggct tctgtctcta tttgccatct	120
taggcataca aaaagagggg tcatcctctt ttacgaattc aatgctacta gggatccaa	180
atactggttg ttgatgactg ccaaaatata ggtatctgct ttcaagaggt catctggctc	240
aaattcaaca tccaatgggg aattttcctg ctctcgaaa ccaaaatat tcagattgta	300
tttgccacgg aggtctaatt tacttcagac tttgacctgc ccaagactga ggaatttca	360
tctccacgat agacacattt ttatccaact gaaagacatc aacctatta tgaaaagaat	420
ggtctgtgct agagactgcc ccatttcata ctctggcgag ataaccgagt cagctccaat	480
cttttctagc actttcttag cggctgact tttgacctta gcaataacag tcggtacccc	540
caaaacttta cagtgcataa ccgcaagcac actcgactcc agattttcac ctgtcgcgac	600
tacaacggta tcgcaggat caatcctgc t	631

<210> SEQ ID NO 31
 <211> LENGTH: 3754
 <212> TYPE: DNA
 <213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 31

ccaatatttt ggtcagcata gtgttctttt tcagtggtaa cagcttgcaa tacttgagca	60
gaaatggcag atttatcaag gaaaaagta acgtaaggtc ctggtgcgac aacttttca	120
aaggcttggc tgttcatttt ttacagccagt tcagccgcaa tcatttggg tgctttacgt	180
tcgacttttg caagagaaaa agcagggaaa gcaatgtctc ccatttctga gtttttaggg	240
gtttccagta actttaaaat agcctcttgg tccaggctat caatgatgct agataattcg	300
ctagcaatca attcttttgt attcattaag agctcctttt tggacttttc tactatttta	360
tcacaatttt aaagaagaa gaaaaattt ttgaaatctc ctgttttttt ggtataatat	420
ggttataaat atagtataa atatagtat aaatatgcac gcaagaggat tttatgagaa	480
aaagagatcg tcatcagtta ataaaaaaaa tgattactga ggagaaatta agtacaaaa	540
aagaattca agatcgggtg gaggcgcaca atgtttgtgt gacgcagaca acctgtctc	600
gtgatttgcg cgaatcggc ttgaccaagg tcaagaaaa tgatatggtg tattatgtac	660
tagtaaatga gacagaaaag attgatttgg tggaaatttt gtctcatcat ttagaagggtg	720
ttgcaagagc agagtttacc ttggtgcttc ataccaaatt gggagaagcc tctgttttgg	780
caaatattgt agatgtaaac aaggatgaat ggattttagg aacagttgct ggtgccaata	840

-continued

ccttattggt tatttgcga gatcagcacg ttgccaaact catggaagat cgtttgctag	900
atttgatgaa agataagtaa ggtcttggga gttgctctca agacttattt ttgaaaagga	960
gagacagaaa atggcgatag aaaagctatc acccggcatg caacagtatg tggatattaa	1020
aaagcaatat ccagatgctt ttttgctctt tcggatgggt gatttttatg aattatttta	1080
tgaggatgcg gtcaatgctg cgcagattct ggaaatttcc ttaacgagtc gcaacaagaa	1140
tgccgacaat ccgatcccta tggcgggtgt tccctatcat tctgcccaac agtatacga	1200
tgtcttgatt gagcagggtt ataaggtggc taccgcagag cagatggaag atcctaaaca	1260
agcagttggg gttgttaaac gagaggttgt tcaggtcatt acgccaggga cagtggtcga	1320
tagcagtaag ccggacagtc agaataattt tttggtttcc atagaccgcg aaggcaatca	1380
atttggccta gcttatatgg atttggtgac gggtgacttt tatgtgacag gtcttttggg	1440
ttcacgctg gtttggggg aaatccgtaa cctcaaggct cgagaagtgg tgttgggtta	1500
tgacttgtct gaggaagaag aacaaatcct cagccgccag atgaatctgg tactctctta	1560
tgaaaaagaa agctttgaa accttcatth attggatttg cgattggcaa cggtgagca	1620
aacggcatct agtaagctgc tccagtatgt tcatcggact cagatgaggg aattgaacca	1680
cctcaaacct gttatccgct acgaaattaa ggatttcttg cagatggatt atgcgaccaa	1740
ggctagtctg gatttgggtg agaatgctcg ctccagtaag aaacaaggca gtcttttctg	1800
gcttttggat gaaacaaaa cggctatggg gatgcgtctc ttgcgttctt ggattcatcg	1860
cccttgatt gataaggaac gaatcgtcca acgtcaagaa gtatgacagg tctttctcga	1920
ccatttcttt gagcgtagtg acttgacaga cagtctcaag ggtgtttatg acattgagcg	1980
cttggctagt cgtgttctt ttggcaaac caatccaaag gatctcttg agttggcgac	2040
tacctgtct agtgtccac ggattcgtgc gattttgaa gggatggagc aacctactct	2100
agcctatctc atcgcaaac tggatgcaat ccctgagttg gagagttga ttagcgcagc	2160
gattgctcct gaagctcctc atgtgattac agatggggga attatccgga ctggatttga	2220
tgagacttta gacaagtatc gttgcgttct cagagaaggg actagctgga ttgctgagat	2280
tgaggctaag gagcgagaaa actctggtat cagcacgctc aagattgact acaataaaaa	2340
ggatgctac tattttcatg tgaccaatc gcaactggga aatgtgccag cccacttttt	2400
ccgcaaggcg acgctgaaaa actcagaacg ctttggaaac gaagaattag cccgtatcga	2460
gggagatag cttgaggcgc gtgagaagtc agccaacctc gaatacgaat tatttatgcg	2520
cattcgtgaa gaggtcggca agtacatcca gcgtttacaa gctctagccc aaggaattgc	2580
gacggttgat gtcttacaga gtctggcggg tgtggctgaa acccagcatt tgattcgacc	2640
tgagtttggg gacgattcac aaattgatat ccgaaaaggg cgcctatgctg tcgttgaaaa	2700
ggttatgggg gctcagacct atattccaaa tacgattcag atggcagaag ataccagtat	2760
tcaattgggt acagggccaa acatgagtg gaaagtctacc tatatgctc agttagccat	2820
gacggcgggt atggcccagc tgggttccta tgttctctgt gaaagcggcc atttaccgat	2880
ttttgatgcg atttttaccg gtatcggagc agcagatgac ttggtttcgg gtcagtcaac	2940
ctttatgggt gagatgatgg aggccaataa tgccatttct catgcgacca agaactctct	3000
cattctcttt gatgaattgg gacgtggaac tgcaacttat gacgggatgg ctcttgctca	3060
gtccatcatc gaatatatcc atgagcacat cggagctaag accctctttg cgaccacta	3120

-continued

```

ccatgagttg actagtctgg agtctagttt acaacacttg gtcaatgtcc acgtggcaac 3180
tttggagcag gatgggcagg tcaccttctc tcacaagatt gaaccgggac cagctgataa 3240
atcctacggt atccatgttg ccaagattgc tggcttgcca gcagaccttt tagcaagggc 3300
ggataagatt ttgactcagc tagagaatca aggaacagag agtcctctc ccatgagaca 3360
aactagtctc gtcactgaac agatttctc ctttgatagg gcagaagagc atcctatcct 3420
agcagaatta gctaaactgg atgtgtataa tatgacacct atgcaggtta tgaatgtctt 3480
agtagagtta aaacagaaa tataaaacca agactcacta gttaatctag ctgtatcaag 3540
gagacttctt tgacaattct ccactttttt gctagaataa catcacacaa acagaatgaa 3600
aagggtgac gcattgtcgc tcccttttgt ctatttttta aggagaaagt atgctgattc 3660
agaaaataaa aacctacaag tggcaggccc tgcttcgctc ctgatgacag gcttgatggt 3720
tgctagtcca cttctgcaac cgcgttatct gcag 3754

```

<210> SEQ ID NO 32

<211> LENGTH: 1337

<212> TYPE: DNA

<213> ORGANISM: Streptococcus pyogenes

<400> SEQUENCE: 32

```

aacaaaaata aagaacttac ctattttcca tccaaaatgt ttagcaatca tcatctgcaa 60
ggcaacgtat tgcattgcat tgatgtgatg agcaactaat atgtcattag aacgttgcgt 120
caaaactagca tctaaataaa gatcgaatg cagttatcaa aaatgcaagc tcctatcggc 180
ccttgtttta attattactc acattgcctt aatgtattta cttgcttatt attactttt 240
ttgctaagtt agtagcgtca gttattcatt gaaaggacat tattatgaaa attcctgtaa 300
caggctttga tccctttggc ggcgaagcta ttaatcctgc ccttgaagct atcaagaaat 360
tgccagcaac cattcatgga gcagaaatca aatgtattga agttccaacg gtttttcaa 420
aatctgccga tgtgctccag cagcatatcg aaagctttca acctgatgca gtcctttgta 480
ttgggcaagc tggtgcccg actggactaa cgccagaacg cgttgccatt aatcaagacg 540
atgctcgcac tcctgataac gaagggaatc agcctattga tacacctatt cgtgcagatg 600
gtaaagcagc ttatttttca accttgccaa tcaaagcgat ggttgctgcc attcatcagg 660
ctgggcttcc tgcttctggt tctaatacag ctggtacett tgtttgcaat cttttgatgt 720
atcaagccct ttacttagtg gataaatatt gtccaaatgc caaagctggg tttatgcata 780
ttccctttat gatggaacag gttgttgata aacctaatc agctgccatg aacctcagat 840
atattacaag aggaattgag gctgctatct ttgccattgt cgatttcaa gatcgttccg 900
atttaaaacg tgtagggggc gctactcact gactgtgacg ctactaaacc tattttaaaa 960
aaacagagat atgaactaac tctgtttttt ttgtgctaaa aatgaaagac ctagggaaac 1020
ttttcatcgg tcttttctca ttgtcatctt aatctaatac tacttctaac atcagcgggt 1080
atagtttgcc agtaattaag aaacgttgtt gatctaaatg agcaatccca ttcaaaacat 1140
taaggtcagc gtaatgggac ttatcaagat ttaaggcttt taacaaagga ctaatatcat 1200
agggtgctac cacctttcca gaatcagggt ggagtttgac aatagtattg gtttgccaaa 1260
tattggcata gagataacca tctacatact ctaattcgtt aagcattgag atagggacac 1320
ttctatagc aactagt 1337

```

-continued

```

<210> SEQ ID NO 33
<211> LENGTH: 1837
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pyogenes

<400> SEQUENCE: 33
tcatgtttga cagcttatca tcgataagct tacttttcga atcaggtcta tccttgaaac    60
aggtgcaaca tagattaggg catggagatt taccagacaa ctatgaacgt atatactcac    120
atcacgcaat cggcaattga tgacattgga actaaattca atcaatttgt tactaacaag    180
caactagatt gacaactaat tctcaacaaa cgttaattta acaacattca agtaactccc    240
accagctcca tcaatgctta cggtaagtaa tcataactta ctaaacctt gttacatcaa    300
ggttttttct ttttgccttg ttcattgagt accataactt tctatattat tgacaactaa    360
attgacaact cttcaattat ttttctgtct actcaaagtt ttcttcattt gatatagtct    420
aattccacca tcactctctc cactctctct accgtcacia cttcatcctc tctcactttt    480
tcgtgtggta acacataatc aaatatcttt cggtttttac gcaactatcg tactgtgtca    540
cctaaaaatat accccttctc aatcgcttct ttaaactcat ctatatataa catatttcat    600
cctcctacct atctattctg aaaaagataa aataaactat tgtttttttt gttattttat    660
aataaaaatta ttaataataag ttaatgtttt ttaaaaatat acaattttat tctatttata    720
gtagctattt ttttcattgt tagtaatatt ggtgaattgt aataaccttt ttaaatctag    780
aggagaacc cagatataaaa tggaggaata ttaatggaaa acaataaaaa agtattgaag    840
aaaaatggtat tttttgtttt agtgacattt cttggactaa caatctcgca agaggatttt    900
gctcaacaag accccgatcc aagccaactt cacagatcta gtttagttaa aaacctcaa    960
aatatatatt ttctttatga gggtgacctt gttactcagc agaatgtgaa atctgttgat    1020
caacttttat ctcacgattt aatatataat gtttcagggc caaattatga taaattaaaa    1080
actgaactta agaaccaaga gatggcaact ttatttaagg ataaaaacgt tgatatttat    1140
ggtgtagaat attaccatct ctgttattta tgtgaaaatg cagaaggagc tgcatgtatc    1200
tacggagggg taacaaatca tgaagggaat catttagaaa ttcctaaaaa gatagtcggt    1260
aaagatcaa tcgatgggat ccaaagccta tcatttgata ttgaacaaa taaaaaatg    1320
gtaactgctc aagaattaga ctataaagt agaaaatctc ttacagataa taagcaacta    1380
tataactaat gaccttctaa atatgaaact ggatatataa agttcatacc taagaataaa    1440
gaaagttttt ggtttgattt tttccctgaa ccagaattta ctcaactcaa atatcttatg    1500
atatataaag ataataaagc gcttgactca aacacaagcc aaattgaagt ctacctaaac    1560
accaagtaac tttttgcttt tggcaacctt acctactgct ggatttagaa attttattgc    1620
aattctttta ttaatgtaaa aaccgctcat ttgatgagcg gttttgtctt atctaaagga    1680
gctttacctc ctaatgctgc aaaattttaa atgttgatt tttgtatttg tctattgtat    1740
ttgatgggta atccatttt tcgacagaca tcgtcgtgcc acctctaaca ccaaaatcat    1800
agacaggagc ttgtagctta gcaactattt tatcgtc                                1837

```

```

<210> SEQ ID NO 34
<211> LENGTH: 841
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pneumoniae

```

-continued

<400> SEQUENCE: 34

```

gatcaatatg tccaagaaac cacatgttcc taagacaaga gctaacagac tggccgctcaa    60
taatagtatt gttctttttt tcatcattac tccttaacta gtgtttaact gattaattag    120
ccagtaaata gtttatcttt atttacacta tctgttaaga tatagtaaaa tgaataaaga    180
acaggacagt caaatcgatt tctaacaatg ttttagaagt agaggatatac tattctaatt    240
tcaatctact atatthttgca cattttcata aaaaaaatga gaactagaac tcacattctg    300
ctctcatttt tcgthtttccc gttctoctat cctgthttta ggagttagaa aatgctgcta    360
cctttactta ctctccttta ataaagccaa tagthtttca gcttctgcca taatagtatt    420
gttgctctgg gtgccaaata gtaaattatt ttttaactct gtgagagtct ctttgccatt    480
ggacttgata attggattct ggatthttcc aagtaaatct tcagcctctc tcagthttct    540
taacctttca gtctcgacct gaggttcttc tgattcctct ggtgattctt ctggtgattc    600
ttcttctggt tcctctgttg gthttggaga ctctggthtc tcgctttgcg gthttctctc    660
tcgagggggt tcttctcag gthtttctgt ctgaggtthc tcctcgthtg gthtttccgt    720
ttgattggta tcagcttgac cattthttgt tctttgaaca tggctgctag cgttaccaaa    780
accattatct gaatgcgacg ttcgthttga tgttcgacat agtacttgac agtcgcaaaa    840
a                                                                 841

```

<210> SEQ ID NO 35

<211> LENGTH: 4500

<212> TYPE: DNA

<213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 35

```

gatcaggaca gtcaaatcga tttctaaca tgthtttagaa gtagatgtgt actattctag    60
tttcaatcta ttatatthtt agaathtttt gthgctagat thgtcaaatt gcttaaaata    120
atthttttca gaaagcaaaa gccgatacct atcgagtagg gtagthcttg ctatcgtcag    180
gcttgctgtg aggtgttaac actthttcaa aatctcttca aacaacgtca gctthgcctt    240
gccgtatata tgttactgac thcgtcagth ctatctgcca cctcaaaacg gtgthttgag    300
ctgacttctg cagthtctatc cacaacctca aaacagtgth ttgagctgac thcgtcagth    360
ctatccacaa cctcaaaaca gtgthttgag ctgactthgt cagthtctatc tacaacctca    420
aaacagtgth ttgagcatca tgcggctagc thcttagtht gctctthgat thtcaattgag    480
tataaaaaca gatgagthtc ththttcttt ttatggacta taaatgttca gctgaaacta    540
ctthcaagga cattattata taaaagaatt thttgaaact aaaatctact atattacact    600
atattgaaag cgtthtthaaa atgaggtata ataaatttac taacacttat aaaaagtgat    660
agaatctatc thtatgtata thtaagata gattgctgta aaaatagtag tagctatgag    720
aaataacaga tagagagaag ggattgaagc ttgaaaaagg ggaataatat gatatttaag    780
gcattcaaga caaaaagca gagaaaaaga caagttgaac tactthttgac agthttthttc    840
gacagthttc tgattgattt atthcttccac thattthgga thgtcccctt taagctggat    900
aagattctga thgtgagctt gattatthtt cccattthtt ctacaagtat thtgccttat    960
gaaaagctat thgaaaaagt gthcagataag gattgagcag gaagtatggt gtaaatagca    1020
taagctgatg tccatcattt gcttataaag agatthttta thtttaattgc agcggthtcc    1080

```


-continued

tggtagataa	actagattgg	caggagtctg	attggagaaa	ggagagggga	aatttggcac	1140
caatttgaga	tagtttgttt	agttcatttt	tgtcatttaa	atgaactgta	gtaaaagaaa	1200
gttaataaaa	gacaaactaa	gtgcattttc	tggaataaat	gtcttatttc	agaaatcggg	1260
atatagatat	agagaggaac	agtatgaatc	ggagtgttca	agaacgtaag	tgctgttata	1320
gcattaggaa	actatcggta	ggagcggttt	ctatgattgt	aggagcagtg	gtatttggaa	1380
cgtctcctgt	tttagctcaa	gaaggggcaa	gtgagcaacc	tctggcaaat	gaaactcaac	1440
tttcggggga	gagctcaacc	ctaactgata	cagaaaagag	ccagccttct	tcagagactg	1500
aactttctgg	caataagcaa	gaacaagaaa	ggaaagataa	gcaagaagaa	aaaattccaa	1560
gagattacta	tgcacgagat	ttggaaaatg	tcgaaacagt	gatagaaaaa	gaagatggtg	1620
aaaccaatgc	ttcaaatggt	cagagagttg	atztatcaag	tgaactagat	aaactaaaga	1680
aacttgaaaa	cgcaacagtt	cacatggagt	ttaagccaga	tgccaaggcc	ccagcattct	1740
ataatctcct	ttctgtgtca	agtgtacta	aaaaagatga	gtacttcact	atggcagttt	1800
acaataatac	tgctactcta	gaggggctg	gttcggatgg	gaaacagttt	tacaataatt	1860
acaacgatgc	acccttaaaa	gttaaacag	gtcagtggaa	ttctgtgact	ttcacagttg	1920
aaaaaccgac	agcagaacta	cctaaaggcc	gagtgcgctt	ctacgtaaac	ggggatttat	1980
ctcgaacaag	tctgagatct	ggcaatttca	ttaaagatat	gccagatgta	acgcatgtgc	2040
aaatcggagc	aaccaagcgt	gccaacaata	cggtttgggg	gtcaaatcta	cagattcggga	2100
atctcactgt	gtataatcgt	gctttaacac	cagaagaggt	acaaaaacgt	agtcaacttt	2160
ttaaagcctc	agatttagaa	aaaaaactac	ctgaaggagc	ggctttaaca	gagaaaacgg	2220
acatatcga	aagcggcgtg	aacggtaaac	caaataaaga	tggaatcaag	agttatcgtg	2280
ttccagcact	tctcaagaca	gataaaggaa	ctttgatcgc	aggtgcagat	gaacgcccgc	2340
tccattcag	tgactggggt	gatatcggta	tggtcatcag	acgtagtcaa	gataatggta	2400
aaacttgggg	tgaccgagta	accattacca	acttacgtga	caatccaaaa	gcttctgacc	2460
catcgatcgg	ttcaccagtg	aatatcgata	tggtgttggg	tcaagatcct	gaaaccaaac	2520
gaatcttttc	tatctatgac	atgttccag	aaggggaagg	aatctttgga	atgtcttcac	2580
aaaaagaaga	agcctacaaa	aaaatcgatg	gaaaaaccta	tcaaatcctc	tatcgtgaag	2640
gagaaaagg	agcttatacc	attcgagaaa	atggtactgt	ctatacacca	gatggtaagg	2700
cgacagacta	tcgctgtgtt	gtagatcctg	ttaaaccagc	ctatagcgac	aagggggatc	2760
tatacaagg	taaccaatta	ctaggcaata	tctacttcac	aacaaacaaa	acttctccat	2820
ttagaattgc	caaggatagc	tatctatgga	tgctctacag	tgatgacgac	gggaagacat	2880
ggtcagcgcc	tcaagatatt	actccgatgg	tcaaagccga	ttggatgaaa	ttcttgggtg	2940
taggtcctgg	aacaggaatt	gtacttcgga	atgggcctca	caagggacgg	attttgatac	3000
cggtttatac	gactaataat	gtatctcact	taaatggctc	gcaatcttct	cgtatcatct	3060
attcagatga	tcattgaaaa	acttggcatg	ctggagaagc	ggtcaacgat	aacogtcagg	3120
tagacggtca	aaagatccac	tcttctacga	tgaacaatag	acgtgcgcaa	aatacagaat	3180
caacggtggt	acaactaaac	aatggagatg	ttaaactctt	tatgcgtggg	ttgactggag	3240
atcttcaggt	tgctacaagt	aaagacggag	gagtgacttg	ggagaaggat	atcaaactgt	3300
atccacaggt	taaagatgtc	tatgttcaaa	tgtctgctat	ccatagcatg	cacgaaggaa	3360

-continued

```

aagaatacat catcctcagt aatgcagggtg gaccgaaacg tgaaaatggg atggtccact 3420
tggcacgtgt cgaagaaaat ggtgagttga cttggctcaa acacaatcca attcaaaaag 3480
gagagtttgc ctataattcg ctccaagaat taggaaatgg ggagtatggc atcttgtatg 3540
aacatactga aaaaggacaa aatgcctata ccctatcatt tagaaaattht aattgggact 3600
ttttgagcaa agatctgatt tctcctaccg aagcgaagt gaagcgaact agagagatgg 3660
gcaaaggagt tattggcttg gagttcgact cagaagtatt ggtcaacaag gctccaaccc 3720
ttcaattggc aaatggtaaa acagcacgct tcatgaccca gtatgatata aaaaccctcc 3780
tatttacagt ggattcagag gatatgggtc aaaaagttac aggtttggca gaaggtgcaa 3840
ttgaaagtat gcataattta ccagtctctg tggcggggcac taagctttcg aatggaatga 3900
acggaagtga agctgctgtt catgaagtgc cagaatacac aggccatta gggacatccg 3960
gcgaagagcc agctccaaca gtcgagaagc cagaatacac aggccacta gggacatccg 4020
gcgaagagcc agccccgaca gtcgagaagc cagaatacac aggccacta gggacagctg 4080
gtgaagaagc agctccaaca gtcgagaagc cagaatttac agggggagtt aatggtacag 4140
agccagctgt tcatgaaatc gcagagtata agggatctga ttcgcttgta actcttacta 4200
caaaagaaga ttatacttac aaagctcctc ttgctcagca ggcacttccct gaaacaggaa 4260
acaaggagag tgacctccta gcttactag gactaacagc tttcttcctt ggtctgttta 4320
cgctagggaa aaagagagaa caataagaga agaattctaa acatttgatt ttgtaaaaat 4380
agaaggagat agcaggtttt caagcctgct atcttttttt gatgacattc aggctgatac 4440
gaaatcataa gaggtctgaa actactttca gagtagtctg ttctataaaa tatagtagat 4500

```

```

<210> SEQ ID NO 36
<211> LENGTH: 705
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus epidermidis

```

```

<400> SEQUENCE: 36

```

```

gatccaagct tatcgatec atcaaaaagt tggcgaacct tttcaaattt tggttcaaat 60
tcttgagatg tatagaatcc aaaatattta ccatttgcat agtctgattg ctcaaagtct 120
tgatactttt ctccacgctc ttttgcaatt tccattgaac gttcagatgga ataatagtcc 180
ataatcataa agaataatatt agcaaagtct tttgcttctt cagattcata gccaaatttt 240
tttttagcta gataaccatg taagttcatt actcctagtc caacagaatg tagttcacta 300
ttcgcttttt ttacacctgg tgcattttga atatttgctt catcacttac aactgtaaga 360
gcatccatac ctgtgaacac agaatctctg aatttacctg attccataac attcactata 420
ttcaatgagc ctaagttaca tgaatatctt cttttaattt catcttcaat tccatagtcg 480
ttaattactg atgtctcttg taattggaaa atttcagtac ataaattact cattttaatt 540
tgccaatat ttgaattcgc atgtactttg tttgcattat ctttaaacat aagatatgga 600
taaccagact gtaattgtgt ttgtgcaatc atatttaaca tttcacgtgc gtcttttttc 660
tttttatcga tttcgaaccc ggggtaccga attcctcgag tctag 705

```

```

<210> SEQ ID NO 37
<211> LENGTH: 442
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus aureus

```

-continued

```

<400> SEQUENCE: 37
gatcaatcctt tgtcgggtaca cgatattcctt cacgactaaa taaacgctca ttcgcgattt    60
tataaatgaa tgttgataac aatgttgat tatctactga aatctcatta cgttgcatcg    120
gaaacattgt gttctgtatg taaaagccgt cttgataatc tttagtagta ccgaagctgg    180
tcatacgaga gttatatttt ccagccaaaa cgatattttt ataatcatta cgtgaaaaag    240
gtttcccttc attatcacac aaatatttta gcttttcagt ttctatatca actgtagctt    300
ctttatccat acgttgaata attgtacgat tctgacgcac catcttttgc acacctttaa    360
tgttatttgt tttaaaagca tgaataagtt tttcaacaca acgatgtgaa tcttctaaga    420
agtcaccgta aaatgaagga tc                                          442

<210> SEQ ID NO 38
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Enterococcus faecalis

<400> SEQUENCE: 38
gcaatacagg gaaaaatgtc                                          20

<210> SEQ ID NO 39
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Enterococcus faecalis

<400> SEQUENCE: 39
cttcatcaaa caattaactc                                          20

<210> SEQ ID NO 40
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Enterococcus faecalis

<400> SEQUENCE: 40
gaacagaaga agccaaaaaa                                          20

<210> SEQ ID NO 41
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Enterococcus faecalis

<400> SEQUENCE: 41
gcaatcccaa ataatacggg                                          20

<210> SEQ ID NO 42
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 42
gctttccagc gtcattattg                                          19

<210> SEQ ID NO 43
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 43

```

-continued

gatctcgaca aaatggtga 19

<210> SEQ ID NO 44
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 44

cacccgcttg cgtggaagc tgccc 25

<210> SEQ ID NO 45
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 45

cgtttggtga ttccagttcc atccg 25

<210> SEQ ID NO 46
<211> LENGTH: 17
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 46

tcacccgctt gcgtggc 17

<210> SEQ ID NO 47
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 47

ggaactggaa tccacaaac 19

<210> SEQ ID NO 48
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 48

tgaagcactg gccgaaatgc tgcgt 25

<210> SEQ ID NO 49
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 49

gatgtacagg attcgttgaa ggctt 25

<210> SEQ ID NO 50
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 50

tagcgaaggc gtagcagaaa ctaac 25

<210> SEQ ID NO 51
<211> LENGTH: 25
<212> TYPE: DNA

-continued

<213> ORGANISM: Escherichia coli
<400> SEQUENCE: 51
gcaaccgaa ctcaacgccc gattt 25

<210> SEQ ID NO 52
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli
<400> SEQUENCE: 52
atacacaagg gtcgcatctg cggcc 25

<210> SEQ ID NO 53
<211> LENGTH: 26
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli
<400> SEQUENCE: 53
tgcgatatgca ttgcagacct tgtggc 26

<210> SEQ ID NO 54
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli
<400> SEQUENCE: 54
gctttcactg gatatcgccg ttggg 25

<210> SEQ ID NO 55
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli
<400> SEQUENCE: 55
gcaaccgaa ctcaacgcc 19

<210> SEQ ID NO 56
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli
<400> SEQUENCE: 56
gcagatgcga cccttgtgt 19

<210> SEQ ID NO 57
<211> LENGTH: 23
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae
<400> SEQUENCE: 57
gtggtgctgt tcagcgcttt cac 23

<210> SEQ ID NO 58
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae
<400> SEQUENCE: 58
gcgatattca caccctacgc agcca 25

-continued

<210> SEQ ID NO 59
<211> LENGTH: 26
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 59

gtcgaaaatg ccggaagagg tatacg 26

<210> SEQ ID NO 60
<211> LENGTH: 26
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 60

actgagctgc agaccggtaa aactca 26

<210> SEQ ID NO 61
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 61

gacagtcagt tcgtcagcc 19

<210> SEQ ID NO 62
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 62

cgtagggtgt gaatatcgc 19

<210> SEQ ID NO 63
<211> LENGTH: 26
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 63

cgtgatggat attcttaacg aagggc 26

<210> SEQ ID NO 64
<211> LENGTH: 23
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 64

accaaactgt tgagccgcct gga 23

<210> SEQ ID NO 65
<211> LENGTH: 23
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 65

gtgatcgccc ctcatctgct act 23

<210> SEQ ID NO 66
<211> LENGTH: 26
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

-continued

<400> SEQUENCE: 66
cgcccttcgt taagaatc catcac 26

<210> SEQ ID NO 67
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 67
tcgcccctca tctgtact 19

<210> SEQ ID NO 68
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 68
gatcgtgatg gatattctt 19

<210> SEQ ID NO 69
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 69
caggaagatg ctgcaccggt tgttg 25

<210> SEQ ID NO 70
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Proteus mirabilis

<400> SEQUENCE: 70
tggttcactg actttgcgat gtttc 25

<210> SEQ ID NO 71
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Proteus mirabilis

<400> SEQUENCE: 71
tcgaggatgg catgcactag aaaat 25

<210> SEQ ID NO 72
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Proteus mirabilis

<400> SEQUENCE: 72
cgctgattag gtttcgctaa aatcttatta 30

<210> SEQ ID NO 73
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Proteus mirabilis

<400> SEQUENCE: 73
ttgatcctca ttttattaat cacatgacca 30

<210> SEQ ID NO 74

-continued

<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 74
gaaacatcgc aaagtcagt 19

<210> SEQ ID NO 75
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 75
ataaaatgag gatcaagttc 20

<210> SEQ ID NO 76
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 76
ccgcctttag cattaattgg tgtttatagt 30

<210> SEQ ID NO 77
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 77
cctattgcag ataccttaaa tgtcttgggc 30

<210> SEQ ID NO 78
<211> LENGTH: 26
<212> TYPE: DNA
<213> ORGANISM: *Streptococcus pneumoniae*

<400> SEQUENCE: 78
agtaaatga aataagaaca ggacag 26

<210> SEQ ID NO 79
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: *Streptococcus pneumoniae*

<400> SEQUENCE: 79
aaaacaggat aggagaacgg gaaaa 25

<210> SEQ ID NO 80
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 80
ttgagtgatg atttactga ctccc 25

<210> SEQ ID NO 81
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: *Proteus mirabilis*

<400> SEQUENCE: 81

-continued

gtcagacagt gatgctgacg acaca 25

<210> SEQ ID NO 82
<211> LENGTH: 27
<212> TYPE: DNA
<213> ORGANISM: Proteus mirabilis

<400> SEQUENCE: 82

tggttgtcat gctgtttgtg tgaaaaat 27

<210> SEQ ID NO 83
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 83

cgagcgggtg gtgttcac 19

<210> SEQ ID NO 84
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 84

caagtcgtcg tcggaggga 19

<210> SEQ ID NO 85
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 85

tcgctgttca tcaagaccc 19

<210> SEQ ID NO 86
<211> LENGTH: 19
<212> TYPE: DNA
<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 86

cggagaacca gacttcac 19

<210> SEQ ID NO 87
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 87

aatgcggtg tacctcggcg ctggt 25

<210> SEQ ID NO 88
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 88

ggcggagggc cagttgcacc tgcca 25

<210> SEQ ID NO 89
<211> LENGTH: 25
<212> TYPE: DNA

-continued

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 89

agccctgctc ctcggcagcc tctgc 25

<210> SEQ ID NO 90

<211> LENGTH: 25

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 90

tggcttttgc aaccgcgttc aggtt 25

<210> SEQ ID NO 91

<211> LENGTH: 25

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 91

gcgcccgcga gggcatgctt cgatg 25

<210> SEQ ID NO 92

<211> LENGTH: 25

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 92

acctgggcgc caactacaag ttcta 25

<210> SEQ ID NO 93

<211> LENGTH: 25

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 93

ggctacgctg ccgggctgca ggccg 25

<210> SEQ ID NO 94

<211> LENGTH: 25

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 94

ccgatctaca ccatcgagat gggcg 25

<210> SEQ ID NO 95

<211> LENGTH: 25

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 95

gagcgcggct atgtgttcgt cggct 25

<210> SEQ ID NO 96

<211> LENGTH: 29

<212> TYPE: DNA

<213> ORGANISM: *Staphylococcus saprophyticus*

<400> SEQUENCE: 96

cgtttttacc cttacctttt cgtactacc 29

-continued

<210> SEQ ID NO 97
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 97

tcaggcagag gtagtacgaa aaggaaggg 30

<210> SEQ ID NO 98
<211> LENGTH: 26
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 98

cgtttttacc cttacctttt cgtact 26

<210> SEQ ID NO 99
<211> LENGTH: 28
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 99

atcgatcatc acattccatt tgttttta 28

<210> SEQ ID NO 100
<211> LENGTH: 27
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 100

caccaagttt gacacgtgaa gattcat 27

<210> SEQ ID NO 101
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 101

atgagtgaag cggagtcaga ttatgtgcag 30

<210> SEQ ID NO 102
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 102

cgctcattac gtacagtgac aatcg 25

<210> SEQ ID NO 103
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 103

ctggtagct tgactcttaa caatcttgct 30

<210> SEQ ID NO 104
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

-continued

<400> SEQUENCE: 104
gacgcgattg tcaactgtacg taatgagcga 30

<210> SEQ ID NO 105
<211> LENGTH: 28
<212> TYPE: DNA
<213> ORGANISM: Haemophilus influenzae

<400> SEQUENCE: 105
gcgtcagaaa aagtaggcga aatgaaag 28

<210> SEQ ID NO 106
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Haemophilus influenzae

<400> SEQUENCE: 106
agcggctcta tcttgtaatg acaca 25

<210> SEQ ID NO 107
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Haemophilus influenzae

<400> SEQUENCE: 107
gaaacgtgaa ctcccctcta tataa 25

<210> SEQ ID NO 108
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Moraxella catarrhalis

<400> SEQUENCE: 108
gccccaaaac aatgaaacat atggt 25

<210> SEQ ID NO 109
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Moraxella catarrhalis

<400> SEQUENCE: 109
ctgcagattt tggaatcata tcgcc 25

<210> SEQ ID NO 110
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Moraxella catarrhalis

<400> SEQUENCE: 110
tggtttgacc agtatttaac gccat 25

<210> SEQ ID NO 111
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Moraxella catarrhalis

<400> SEQUENCE: 111
caacggcacc tgatgtacct tgtac 25

<210> SEQ ID NO 112

-continued

<211> LENGTH: 18
<212> TYPE: DNA
<213> ORGANISM: *Moraxella catarrhalis*

<400> SEQUENCE: 112
ggcacctgat gtaccttg 18

<210> SEQ ID NO 113
<211> LENGTH: 17
<212> TYPE: DNA
<213> ORGANISM: *Moraxella catarrhalis*

<400> SEQUENCE: 113
aacagctcac acgcatt 17

<210> SEQ ID NO 114
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: *Moraxella catarrhalis*

<400> SEQUENCE: 114
ttacaacctg caccacaagt catca 25

<210> SEQ ID NO 115
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: *Moraxella catarrhalis*

<400> SEQUENCE: 115
gtacaaacaa gccgtcagcg actta 25

<210> SEQ ID NO 116
<211> LENGTH: 23
<212> TYPE: DNA
<213> ORGANISM: *Moraxella catarrhalis*

<400> SEQUENCE: 116
caatctgcgt gtgtgcgcttc act 23

<210> SEQ ID NO 117
<211> LENGTH: 26
<212> TYPE: DNA
<213> ORGANISM: *Moraxella catarrhalis*

<400> SEQUENCE: 117
gctactttgt cagcttttagc cattca 26

<210> SEQ ID NO 118
<211> LENGTH: 24
<212> TYPE: DNA
<213> ORGANISM: *Moraxella catarrhalis*

<400> SEQUENCE: 118
tgttttgagc tttttatttt ttga 24

<210> SEQ ID NO 119
<211> LENGTH: 22
<212> TYPE: DNA
<213> ORGANISM: *Moraxella catarrhalis*

<400> SEQUENCE: 119

-continued

cgctgacggc ttgtttgtac ca 22

<210> SEQ ID NO 120
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 120

tctgtgctag agactgcccc atttc 25

<210> SEQ ID NO 121
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 121

cgatgtcttg attgagcagg gttat 25

<210> SEQ ID NO 122
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 122

atcccacctt aggcggctgg ctcca 25

<210> SEQ ID NO 123
<211> LENGTH: 31
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 123

acgtcaagtc atcatggccc ttacgagtag g 31

<210> SEQ ID NO 124
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 124

gtgtgacggg cggtgtgtac aaggc 25

<210> SEQ ID NO 125
<211> LENGTH: 28
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 125

gagttgcaga ctccaatccg gactacga 28

<210> SEQ ID NO 126

-continued

<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 126

ggaggaaggt ggggatgacg 20

<210> SEQ ID NO 127
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 127

atggtgtgac gggcgggtg 20

<210> SEQ ID NO 128
<211> LENGTH: 32
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 128

ccctatacat caccttgcg ttagcagag ag 32

<210> SEQ ID NO 129
<211> LENGTH: 28
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 129

ggggggacca tcctccaag ctaaatac 28

<210> SEQ ID NO 130
<211> LENGTH: 32
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Description of Artificial Sequence:
Oligonucleotide

<400> SEQUENCE: 130

cgtccacttt cgtgtttgca gaggctgtg tt 32

<210> SEQ ID NO 131
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 131

caggagtacg gtgattttta 20

<210> SEQ ID NO 132
<211> LENGTH: 20

-continued

<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 132
atttctggtt tggtcataca 20

<210> SEQ ID NO 133
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Proteus mirabilis

<400> SEQUENCE: 133
cgggagtcag tgaatcatc 20

<210> SEQ ID NO 134
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Proteus mirabilis

<400> SEQUENCE: 134
ctaaaatcgc cacacctctt 20

<210> SEQ ID NO 135
<211> LENGTH: 18
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 135
gcagcgtggt gtcgttca 18

<210> SEQ ID NO 136
<211> LENGTH: 18
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 136
agctggcaac ggctggtc 18

<210> SEQ ID NO 137
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 137
attcacacc tacgcagcca 20

<210> SEQ ID NO 138
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

<400> SEQUENCE: 138
atccgcagc atctctttgt 20

<210> SEQ ID NO 139
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 139
ctggttagct tgactcttaa caatc 25

-continued

<210> SEQ ID NO 140
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus saprophyticus

<400> SEQUENCE: 140
tcttaacgat agaatggagc aactg 25

<210> SEQ ID NO 141
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pyogenes

<400> SEQUENCE: 141
tgaaaattct tgtaacaggc 20

<210> SEQ ID NO 142
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pyogenes

<400> SEQUENCE: 142
ggccaccagc ttgcccaata 20

<210> SEQ ID NO 143
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pyogenes

<400> SEQUENCE: 143
atattttctt tatgagggtg 20

<210> SEQ ID NO 144
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pyogenes

<400> SEQUENCE: 144
atccttaa at aaagttgcca 20

<210> SEQ ID NO 145
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus epidermidis

<400> SEQUENCE: 145
atcaaaaagt tggcgaacct tttca 25

<210> SEQ ID NO 146
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus epidermidis

<400> SEQUENCE: 146
caaaagagcg tggagaaaag tatca 25

<210> SEQ ID NO 147
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus epidermidis

-continued

<400> SEQUENCE: 147
tctcttttaa tttcatcttc aattccatag 30

<210> SEQ ID NO 148
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus epidermidis

<400> SEQUENCE: 148
aaacacaatt acagtctggt tatccatctc 30

<210> SEQ ID NO 149
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus aureus

<400> SEQUENCE: 149
cttcatttta cggtgacttc ttagaagatt 30

<210> SEQ ID NO 150
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus aureus

<400> SEQUENCE: 150
tcaactgtag cttctttatc catacgttga 30

<210> SEQ ID NO 151
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus aureus

<400> SEQUENCE: 151
atattttagc ttttcagttt ctatatcaac 30

<210> SEQ ID NO 152
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus aureus

<400> SEQUENCE: 152
aatctttgtc ggtacacgat attcttcacg 30

<210> SEQ ID NO 153
<211> LENGTH: 30
<212> TYPE: DNA
<213> ORGANISM: Staphylococcus aureus

<400> SEQUENCE: 153
cgtaatgaga tttcagtaga taatacaaca 30

<210> SEQ ID NO 154
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Haemophilus influenzae

<400> SEQUENCE: 154
tttaacgatc cttttactcc ttttg 25

-continued

```

<210> SEQ ID NO 155
<211> LENGTH: 25
<212> TYPE: DNA
<213> ORGANISM: Haemophilus influenzae

<400> SEQUENCE: 155

actgctgttg taaagaggtt aaaat                25

<210> SEQ ID NO 156
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 156

atttggtgac gggtgacttt                    20

<210> SEQ ID NO 157
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 157

gctgaggatt tgttcttctt                    20

<210> SEQ ID NO 158
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 158

gagcggtttc tatgattgta                    20

<210> SEQ ID NO 159
<211> LENGTH: 20
<212> TYPE: DNA
<213> ORGANISM: Streptococcus pneumoniae

<400> SEQUENCE: 159

atctttcctt tcttgttctt                    20

<210> SEQ ID NO 160
<211> LENGTH: 18
<212> TYPE: DNA
<213> ORGANISM: Moraxella catarrhalis

<400> SEQUENCE: 160

gctcaaatca gggtcagc                      18

<210> SEQ ID NO 161
<211> LENGTH: 861
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 161

atgagtattc aacatttcg tgtcgccctt attccctttt ttgcggcatt ttgccttctt    60
gtttttgctc acccagaaac gctggtgaaa gtaaaagatg ctgaagatca gttgggtgca    120
cgagtggggtt acatcgaact ggatctcaac agcggtaaga tccttgagag ttttcgcccc    180
gaagaacggt ttccaatgat gagcactttt aaagtcttgc tatgtggcgc ggtattatcc    240
cgtgttgacg ccgggcaaga gcaactcggg cgccgcatac actatttctca gaatgacttg    300

```

-continued

```

gttgagtact caccagtcac agaaaagcat cttacggatg gcatgacagt aagagaatta 360
tgcagtgctg ccataacat gagtgataac actgctggcca acttacttct gacaacgatc 420
ggaggaccga aggagctaac cgcttttttg cacaacatgg gggatcatgt aactcgcctt 480
gatcgttggg aaccggagct gaatgaagcc ataccaaaacg acgagcgtga caccacgatg 540
cctgcagcaa tggcaacaac gttgcgcaaa ctattaactg gcgaactact tactctagct 600
tcccggcaac aattaataga ctggatggag gcggataaag ttgcaggacc acttctgcgc 660
tcggcccttc cggctggctg gtttattgct gataaatctg gagccggtga gcgtgggtct 720
cgcggtatca ttgcagcact ggggccagat ggtaagccct cccgtatcgt agttatctac 780
acgacgggga gtcaggcaac tatggatgaa cgaaatagac agatcgctga gataggtgcc 840
tcactgatta agcattggta a 861

```

```

<210> SEQ ID NO 162
<211> LENGTH: 918
<212> TYPE: DNA
<213> ORGANISM: Pasteurella haemolytica

```

```

<400> SEQUENCE: 162

```

```

atgttaaata agttaaaat cggcacatta ttattgctga cattaacggc ttgttcgccc 60
aattctgttc attcggtaac gtctaatacgc cagcctgcta gtgcgcctgt gcaacaatca 120
gccacacaag ccacctttca acagactttg gcgaatttgg aacagcagta tcaagcccga 180
attggcgttt atgtatggga tacagaaaacg ggacattcct tgtcttatcg tgcagatgaa 240
cgctttgctt atgcgtccac tttcaaggcg ttgttggctg gggcgggtgt gcaatcgctg 300
cctgaaaaag atttaaatcg taccatttca tatagccaaa aagatttggg tagttattct 360
cccgaaaccc aaaaatacgt tggcaaaagc atgacgattg cccaattatg tgaagcagcc 420
gtcgggttta gcgacaacag cgcgaccaat ttgctgctca aagaattggg tggcgtggaa 480
caatatcaac gtattttgcg acaattaggc gataacgtaa cccataccaa tcggctagaa 540
cccgatthaa atcaagccaa acccaacgat attcgtgata cgagtacacc caaacaatg 600
gcgatgaatt taaatcgcta tttattgggc aacacattaa ccgaatcgca aaaaacgatt 660
ttgtggaatt gtttgacaa taacgcaaca ggcaatccat tgattcgcgc tgctacgcca 720
acatcgtgga aagtgtacga taaaagcggg gcgggtaaat atggtgtacg caatgatatt 780
gcggtggttc gcatacaaa tcgcaaacgc attgtgatgg caatcatgag tacgcaattt 840
accgaagaag ccaaattcaa caataaatta gtagaagatg cagcaaagca agtatttcat 900
actttacagc tcaactaa 918

```

```

<210> SEQ ID NO 163
<211> LENGTH: 864
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

```

```

<400> SEQUENCE: 163

```

```

atgogttata ttgcctgtg tattatctcc ctgtagcca ccctgccgct ggcggtacac 60
gccagcccgc agccgcttga gcaaattaaa ctaagcgaag gccagctgtc gggccgcgta 120
ggcatgatag aatggatct gccacggcg cgcacgctga ccgctggcg cgcgcatgaa 180
cgctttccca tgatgagcac ctttaaagta gtgctctgcg gcgcagtgtc ggcgcggtg 240

```

-continued

```

gatgccggtg acgaacagct ggagcgaaa atccactatc gccagcagga tctggtggac 300
tactcgccgg tcagcgaaaa acaccttgcc gacgcaatga cggtcggcga actctcgccc 360
gccgccatta ccatgagcga taacagcgcc gccaatctgc tactggccac cgtcgccggc 420
cccgcaggat tgactgcctt ttgcccag atcggcgaca acgtcacccg ccttgaccgc 480
tgggaaacgg aactgaatga ggcgcttccc ggcgacgccc gcgacaccac taccgccggc 540
agcatggccg cgaccctgcg caacgttggc ctgaccagcc agcgtctgag cggccgttcg 600
caacggcagc tgctgcagtg gatggtggac gatcgggtcg ccggaccggt gatccgctcc 660
gtgctgccgg cgggctggtt tatcgccgat aagaccggag ctggcgagcg gggtcgcggc 720
gggattgtcg cctgctttgg cccgaataac aaagcagagc gcattgtggt gatttatctg 780
cgggatcccc cggcgagcat ggccgagcga aatcagcaaa tcgccgggat cggcaaggcg 840
ctgtacgagc actggcaacg ctaa 864

```

```

<210> SEQ ID NO 164
<211> LENGTH: 534
<212> TYPE: DNA
<213> ORGANISM: Klebsiella pneumoniae

```

<400> SEQUENCE: 164

```

atggacacaa cgcaggtcac attgatacac aaaattctag ctgcggcaga tgagcgaat 60
ctgccctctt ggatcgggtg gggctgggcy atcgatgcac ggctagggcg tgtaacacgc 120
aagcacgatg atattgatct gacgtttccc ggcgagagcg gcggcgagct cgaggcaata 180
gttgaatgc tcggcgggcy cgtcatggag gagttggact atggattctt agcggagatc 240
ggggatgagt tacttgactg cgaacctgct tgggtggcag acgaagcgtg tgaatcgcg 300
gaggctccgc agggctcgtg cccagagggc gctgagggcg tcatcgccgg gcggccagtc 360
cgttgtaaca gctgggagc gatcatctgg gattactttt actatgccga tgaagtacca 420
ccagtggact ggctacaaa gcacatagag tctacagggc tcgcatgcac ctactcggg 480
gcggaagg ttgaggtctt gcgtgccgct ttcaggtcgc gatatgcggc ctaa 534

```

```

<210> SEQ ID NO 165
<211> LENGTH: 465
<212> TYPE: DNA
<213> ORGANISM: Unknown Organism
<220> FEATURE:
<223> OTHER INFORMATION: Description of Unknown Organism:
    Enterobacteriaceae

```

<400> SEQUENCE: 165

```

atgggcatca ttcgcacatg taggctcggc cctgaccaag tcaaatccat gcgggctgct 60
cttgatcttt tcggctcgtg gttcggagac gtagccacct actccaaca tcagccggac 120
tccgattacc tcgggaactt gctccgtagt aagacattca tcgcgcttgc tgccttcgac 180
caagaagcgg ttgttgccgc tctcgggct tacgttctgc ccaggttga gcagccgct 240
agtgagatct atatctatga tctcgcagtc tccggcgagc accggaggca gggcattgcc 300
accgcgctca tcaatctcct caagcatgag gccaacgcgc ttggtgctta tgtgatctac 360
gtgcaagcag attacggtga cgatcccgca gtggctctct atacaaagtt gggcatacgg 420
gaagaagtga tgcactttga tatcgaccca agtaccgcca ctaa 465

```

-continued

```

<210> SEQ ID NO 166
<211> LENGTH: 861
<212> TYPE: DNA
<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 166
atgcatacgc ggaagcaat aacggaggcg cttcaaaac tcggagtcca aaccggtgac   60
ctattgatgg tgcatgcctc acttaaagcg attggtccgg tcgaaggagg agcgggagacg   120
gtcgttgccg cgttacgctc cgcggttggg ccgactggca ctgtgatggg atacgcatcg   180
tgggaccgat caccctacga ggagactcgt aatggcgctc ggttgatga caaaaccgcg   240
cgtacctggc cgccgttcga tcccgcacg gccgggactt accgtgggtt cggcctgctg   300
aatcagtttc tggttcaagc ccccggcgcg cggcgagcg cgcaccccg tgcacgatg   360
gtcgcggttg gtccactggc tgaacgctg acggagcctc acaagctcgg tcacgccttg   420
ggggaagggt cgcccgctga gcggttcggt cgccttgccg ggaaggccct gctggtgggt   480
gcccgcgctaa actccgttac cgcattgcac taccccgagg cgggtgcccga tatcccac   540
aaacggcggg tgacgtatga gatgcccgat cttggaagca acggcgaagt cgcctggaaa   600
acggcatcgg attacgattc aaacggcatt ctcgattgct ttgctatcga aggaaagccg   660
gatgcggtcg aaactatagc aaatgcttac gtgaagctcg gtcgccatcg agaaggtgtc   720
gtgggctttg ctcaagtcta cctgttcgac gcgcaggaca tcgtgacggt cggcgtcacc   780
tatcttgaga agcatttcgg aaccactcgg atcgtgccag cacacgaagt cgcgagtg   840
tcttgcgagc cttcaggtta g                                     861

```

```

<210> SEQ ID NO 167
<211> LENGTH: 816
<212> TYPE: DNA
<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 167
atgaccgatt tgaatatccc gcatacacac gcgcacctg tagacgcatt tcaggcgctc   60
ggcatccgcg cggggcagcg gctcatgctg cacgcatccg ttaaagcagt gggcgcggtg   120
atgggcggcc ccaatgtgat cttgcaggcg ctcatggatg cgctcacgcc cgacggcaccg   180
ctgatgatgt atgcccgatg gcaagacatc cccgacttta tcgactcgtt gccggacgcg   240
ctcaaggccg tgtatcttga gcagcaccga ccttttgacc ccgccaccgc ccgcccgtg   300
cgcgaaaaca gcgtgctagc ggaatttttg cgcacatggc cgtgctgca tcgcagcgca   360
aaccocgaag cctctatggt ggcggtagc aggcaggccg ctttgcctgac cgctaatac   420
gcgctggatt atggctacgg agtcgagtcg ccgctggcta aactggtggc aatagaagga   480
tacgtgctga tgcttgccgc gccgctgat accatcacac tgctgcacca cgcggaatat   540
ctggccaaga tgcgccaaa gaacgtggtc cgctaccctg gcccgattct gcgggacggg   600
cgcaaagtgt gggtagccgt tgaggactat gacaccggtg atccgcaaga cgattatagt   660
tttgagcaaa tcgcgcgcga ttatgtggcg caggcggcg gcacacgcg caaagtcggt   720
gatgcccgat cttacctggt cgccgcgag gacctcacac ggtttgcggt gcagtggtt   780
gaatcacggt tcggtgactc agcgtcatic ggatag                                     816

```

```

<210> SEQ ID NO 168

```

-continued

<211> LENGTH: 498

<212> TYPE: DNA

<213> ORGANISM: *Pseudomonas aeruginosa*

<400> SEQUENCE: 168

```

atgctctatg agtggctaaa tcgatctcat atcgtcgagt ggtggggcgg agaagaagca    60
cgcccgacac ttgctgacgt acaggaacag tacttgccaa gcgttttagc gcaagagtcc    120
gtcactccat acattgcaat gctgaatgga gagccgattg ggtatgccca gtcgtacggt    180
gctcttgtaa gcggggacgg atgggtggaa gaagaaaccg atccaggagt acgcggaata    240
gaccagttac tggcgaatgc atcacaactg ggcaaaggct tgggaaccaa gctgggtcga    300
gctctgggtg agttgctgtt caatgatccc gaggtcacca agatccaaac ggaccctgctg    360
ccgagcaact tgcgagcgat ccgatgctac gagaaagcgg ggtttgagag gcaaggtacc    420
gtaaccaccc cagatggtcc agccgtgtac atggttcaaa cacgccaggc attcagcgca    480
acacgcagtg atgcctaa                                498

```

<210> SEQ ID NO 169

<211> LENGTH: 2007

<212> TYPE: DNA

<213> ORGANISM: *Staphylococcus aureus*

<400> SEQUENCE: 169

```

atgaaaaaga taaaaattgt tccacttatt ttaatagttg tagttgtcgg gtttggata    60
tatttttatg cttcaaaaga taaagaaatt aataatacta ttgatgcaat tgaagataaa    120
aatttcaaac aagtttataa agatagcagt tatatttcta aaagcgataa tggatgaagta    180
gaaatgactg aacgtccgat aaaaatatat aatagtttag gcgttaaaga tataaacatt    240
caggatcgta aaataaaaaa agtatctaaa aataaaaaac gagtagatgc tcaatataaa    300
attaaaacaa actacggtaa cattgatcgc aacgttcaat ttaattttgt taaagaagat    360
ggtatgtgga agttagattg ggatcatagc gtcattattc caggaatgca gaaagaccaa    420
agcatacata ttgaaaaattt aaaatcagaa cgtggtaaaa ttttagaccg aaacaatgtg    480
gaattggcca atacaggaac acatatgaga ttaggcacgc ttccaagaa tgtatctaaa    540
aaagattata aagcaatcgc taaagaacta agtatttctg aagactatat caacaacaaa    600
tggatcaaaa ttgggtcaaa gatgatacct tcgttccact ttaaaaccgt taaaaaatg    660
gatgaatatt taagtgattt cgcaaaaaaa tttcatctta caactaatga aacagaaagt    720
cgtaactatc ctctagaaaa agcgacttca catctattag gttatgttgg tcccattaac    780
tctgaagaat taaaacaaaa agaataataa ggctataaag atgatgcagt tattggtaaa    840
aagggactcg aaaaacttta cgataaaaaa ctccaacatg aagatggcta tcgtgtcaca    900
atcgttgacg ataatagcaa tacaatcgca catacattaa tagagaaaa gaaaaaagat    960
ggcaagata ttcaactaac tattgatgct aaagttcaaa agagtattta taacaacatg    1020
aaaaatgatt atggctcagg tactgctatc caccctcaaa caggtgaatt attagcactt    1080
gtaagcacac cttcatatga cgtctatcca tttatgtatg gcatgagtaa cgaagaatat    1140
aataaattaa ccgaagataa aaaagaacct ctgctcaaca agttccagat tacaacttca    1200
ccaggttcaa ctcaaaaaat attaacagca atgattgggt taaataacaa aacattagac    1260
gataaaacaa gttataaaat cgatggtaaa ggttggcaaa aagataaatc ttggggtggt    1320

```

-continued

tacaacgtta caagatatga agtggtaaat ggtaatatcg acttaaaaca agcaatagaa	1380
tcacagata acattttctt tgctagagta gcaactcgaat taggcagtaa gaaattgaa	1440
aaaggcatga aaaaactagg tgttggtaaa gataaccaa gtgattatcc attttataat	1500
gctcaaattt caaacaaaaa tttagataat gaaatattat tagctgattc aggttacgga	1560
caaggtgaaa tactgattaa cccagtacag atcctttcaa tctatagcgc attagaaaaat	1620
aatggcaata ttaacgcacc tcacttatta aaagacacga aaaacaaagt ttggaagaaa	1680
aatattattt ccaaagaaaa tatcaatcta ttaaagtatg gtatgcaaca agtcgtaaat	1740
aaaacacata aagaagatat ttatagatct tatgcaaact taattggcaa atccggtact	1800
gcagaactca aatgaaaca aggagaaagt ggacagacaaa ttgggtgggt tatatcatat	1860
gataaagata atccaaacat gatgatggct attaagtta aagatgtaca agataaagga	1920
atggctagct acaatgccaa aatctcaggt aaagtgtatg atgagctata tgagaacggt	1980
aataaaaaat acgatataga tgaataa	2007

<210> SEQ ID NO 170

<211> LENGTH: 2607

<212> TYPE: DNA

<213> ORGANISM: Enterococcus faecium

<400> SEQUENCE: 170

atgaataaca tcggcattac tgtttatgga tgtgagcagg atgaggcaga tgcattccat	60
gctctttcgc ctgcgctttg cgttatggca acgataatta acgccaacgt gtcggaatcc	120
aacgccaaat ccgcgccttt caatcaatgt atcagtgtgg gacataaatc agagatttcc	180
gcctctattc ttcttgcgct gaagagagcc ggtgtgaaat atatttctac ccgaagcatc	240
ggctgcaatc atatagatag aactgctgct aagagaatgg gcatcactgt cgacaatgtg	300
gcgtactcgc cggatagcgt tgccgattat actatgatgc taattcttat ggcagtacgc	360
aacgtaaaat cgattgtcgc ctctgtgtaa aaacatgatt tcaggttgga cagcgaccgt	420
ggcaaggtag tcagcgacat gacagttggt gtggtgggaa cgggccagat aggcaaagcg	480
gttattgagc ggctgcgagg atttggatgt aaagtgttgg cttatagtcg cagccgaagt	540
atagaggtaa actatgtacc gtttgatgag ttgctgcaaa atagcgatat cgttacgctt	600
catgtgccgc tcaatcgcga tacgcactat attatcagcc acgaacaaat acagagaatg	660
aagcaaggag catttcttat caatactggg cgcggtcac ttgtagatag ctatgagttg	720
gttaaagcat tagaaaacgg gaaactgggc ggtgccgcat tggatgtatt ggaaggagag	780
gaagagtttt tctactctga ttgcacccaa aaaccaattg ataatcaatt tttacttaa	840
cttcaaaaga tgctaactg gataatcaca cgcatacgg octattatac cgagcaagcg	900
ttgcgtgata ccggtgaaaa aaccatataa aactgtttgg attttgaaag gagacaggag	960
catgaataga ataaaagttg caatactggt tgggggttgc tcagaggagc atgacgtatc	1020
ggtaaaatct gcaatagaga tagccgctaa cattaataaa gaaaaatagc agcogttata	1080
cattggaatt acgaaatctg gtgtatgtaa aatgtgcgaa aaaccttgcg cggaatggga	1140
aaacgacaat tgctattcag ctgtactctc gccggataaa aaaatgcacg gattacttgt	1200
taaaaagaac catgaatag aaatcaacca tgttgatgta gcattttcag ctttgcattg	1260
caagtcaggt gaagatggat ccatacaagg tctgtttgaa ttgtccggtg tcccttttgt	1320

-continued

```

aggctgcgat attcaaaagt cagcaatttg tatggacaaa tcgttgacat acatcgttgc 1380
gaaaaatgct gggatagcta ctcccgcott ttgggttatt aataaagatg ataggccggt 1440
ggcagctacg tttacctatc ctggttttgt taagccggcg cgttcaggct catccttcgg 1500
tgtgaaaaaa gtcaaatagcg cggacgaatt ggactacgca attgaatcgg caagacaata 1560
tgacagcaaa atcttaattg agcaggctgt ttcgggctgt gaggtcggtt gtgcggtatt 1620
gggaaacagt gccgcgtagg ttggtggcga ggtggaccaa atcaggctgc agtacggaat 1680
ctttcgtatt catcaggaag tcgagccgga aaaaggctct gaaaacgcag ttataaccgt 1740
tcccgcagac ctttcagcag aggagcggag acggatacag gaaacggcaa aaaaaatata 1800
taaagcgctc ggctgtagag gtctagcccg tgtggatag tttttacaag ataacggccg 1860
cattgtactg aacgaagtca atactctgcc cggtttcacg tcatacagtc gttatccccg 1920
tatgatggcc gctgcaggta ttgcacttcc cgaactgatt gaccgctga tcgtattagc 1980
gttaaagggg tgataagcat ggaaatagga tttacttttt tagatgaaat agtacacggt 2040
gttcgttggg acgctaaata tgccacttgg gataatttca ccgaaaacc ggttgacggt 2100
tatgaagtaa atcgattgt agggacatac gagttggctg aatcgctttt gaaggcaaaa 2160
gaactggctg ctaccaagg gtacggattg cttctatggg acggttaccg tcctaagcgt 2220
gctgtaaaat gttttatgca atgggctgca cagccggaaa ataacctgac aaaggaaagt 2280
tattatccca atattgaccg aactgagatg atttcaaaag gatagctggc ttcaaaatca 2340
agccatagcc gcggcagtg cattgatctt acgctttatc gattagacac gggtagctt 2400
gtaccaatgg ggagccgatt tgattttatg gatgaacgct ctcacatgac ggcaaatgga 2460
atatcatgca atgaagcgca aaatcgaca cgtttgcgct ccatcatgga aaacagtggg 2520
tttgaagcat atagcctcga atggtggcac tatgtattaa gagacgaacc atacccaat 2580
agctattttg atttccccgt taaataa 2607

```

<210> SEQ ID NO 171

<211> LENGTH: 1288

<212> TYPE: DNA

<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 171

```

ggatccatca ggcaacgacg ggctgctgcc ggccatcagc ggacgcaggg aggactttcc 60
gcaaccggcc gttcgatgcg gcaccgatgg ccttcgcgca gggtagtga atccgccagg 120
attgacttgc gctgcacctac ctctcactag tgaggggicgg cagcgcatca agcggtgagc 180
gcactccggc accgccaact ttcagcacat gcgtgtaaat catcgtcgtg gagacgtcgg 240
aatggccgag cagatcctgc acggttcgaa tgcgtaacc gctgcggagc aaggccgtcg 300
cgaacgagtg gcggaggggtg tgcgggtgtg cgggcttcgt gatgcctgct tgttctacgg 360
cacgtttgaa ggcgcgctga aaggctcgtg catacatggt atggcgacgc acgacaccgc 420
tccgtggatc ggtcgaatgc gtgtgctgcg caaaaaccca gaaccacggc caggaatgcc 480
cggcgcgicgg atacttccgc tcaagggcgt cgggaagcgc aacgccgctg cggccctcgg 540
cctggtcctt cagccaccat gcccgtcac gcgacagctg ctcgcgcagg ctgggtgcca 600
agctctcggg taacatcaag gcccgatcct tggagccctt gccctccgc acgatgatcg 660
tgccgtgatc gaaatccaga tccttgacct gcagttgcaa accctcactg atccgatgc 720

```

-continued

```

ccgtccata cagaagctgg gcgaacaaac gatgctcgcc ttccagaaaa ccgaggatgc 780
gaaccacttc atccggggtc agcaccaccg gcaagcgccg cgacggccga ggtcttccga 840
tctctgaag ccagggcaga tccgtgcaca gcaccttgcc gtagaagaac agcaaggccg 900
ccaatgcctg acgatgcgtg gagaccgaaa ccttgcgctc gttcgccagc caggacagaa 960
atgctctgac ttcgctgctg cccaaggttg ccgggtgacg cacaccgtgg aaacggatga 1020
aggcacgaac ccagtggaca taagcctgtt cggttcgtaa gctgtaatgc aagtagcgtg 1080
tgcgctcagc caactggttc agaacctga ccgaacgcag cgggtgtaac ggcgcagtgg 1140
cggtttcat ggcttgttat gactgttttt ttgtacagtc tatgcctcgg gcatccaagc 1200
agcaagcgcg ttacgcccgt ggtcgatggt tgatgttatg gagcagcaac gatgttacgc 1260
agcagggcag tcgccctaaa acaaagtt 1288

```

<210> SEQ ID NO 172

<211> LENGTH: 1650

<212> TYPE: DNA

<213> ORGANISM: Pseudomonas aeruginosa

<400> SEQUENCE: 172

```

gtagatgca ctaagcacat aattgctcac agccaaacta tcaggtaag tctgctttta 60
ttatttttaa gcgtgcataa taagccctac acaaattggg agatataatca tgaaggttg 120
gctttttcct gttatcgcaa tagttggcga agtaatcgca acatccgat taaaatctag 180
cgagggcttt actaagcttg ccccttccgc cgttgcata atcggttatg gcatcgcatt 240
ttattttctt tctctggttc tgaatccat cctgtcgggt gttgcttatg cagtctggtc 300
gggactcggc gtcgctataa ttacagccat tgcctggttg cttcatgggc aaaagcttga 360
tgcgtggggc tttgtaggta tggggctcat aattgctgcc tttttgctcg cccgatcccc 420
atcgtggaag tcgctcgga ggcgcagccc atggtgacgg tgttcggcat tctgaatctc 480
accgaggact ctttcttcca tgagagccgg cggctagacc ccgcccggcg tgtcaccgcg 540
gcgatcgaat tgctgcgagt cggatcagac gtcgtggatg tcggaccggc cgccagccat 600
ccggacgcga ggctgtatc gccggccgat gagatcagac gtattgcgcc gctcttagac 660
gccctgtccg atcagatgca ccgtgtttca atcgacagct tccaaccgga aaccagcgcg 720
tatgcgctca agcgcggcgt gggctacctg aacgatatcc aaggatttcc tgaccctgcg 780
ctctatcccc atattgctga ggcggactgc aggtggtgg ttatgcactc agcgcagcgg 840
gatggcatcg ccaccgcac cggtcacctt cgaccggaag acgctcga cgagattgtg 900
cggttcttgc aggcgcgggt ttccgccttg cgacggagcg ggtcgcctgc cgaccgctc 960
atcctcgatc cggggatggg attttcttg agcccgcac cggaaacatc gctgcacgtg 1020
ctgtcgaacc ttcaaaagct gaagtcggcg ttggggcttc cgctattggt ctcggtgtcg 1080
cggaaatcct tcttggggcg caccgttggc cttcctgtaa aggatctggg tccagcgagc 1140
cttgcggcgg aacttcacgc gatcggcaat ggcgctgact acgtccgcac ccaocgcct 1200
ggagatctgc gaagcgaat caccttctcg gaaaccctcg cgaaatttcg cagtccgcac 1260
gccagagacc gagggtaga tcatgcttag cattcacctt ccggccgccc gctagcggac 1320
cctggtcagg ttccgcgaag gtgggcgag acatgctggg ctcgtcagga tcaaaactgca 1380
ctatgaggcg gcggttcata ccgcgccagg ggagcgaatg gacagcgagg agcctccgaa 1440

```

-continued

```

cgttcggggtc gcctgctcgg gtgatatcga cgagggttggtg cggctgatgc acgacgctgc 1500
ggcgtggatg tccgccaagg gaacgcccg ctagggacgtc gcgctgatgc accggacatt 1560
cgcggagacc ttcgtcctga gatccgagct cctagtcgcg agttgcagcg acggcatcgt 1620
cggctgttgc accttgtcgg ccgaggatcc 1650

```

```

<210> SEQ ID NO 173
<211> LENGTH: 630
<212> TYPE: DNA
<213> ORGANISM: Enterococcus faecium

```

```
<400> SEQUENCE: 173
```

```

atgggtccga atcctatgaa aatgtatcct atagaaggaa acaaatcagt acaatttattc 60
aaacctattt tagaaaaatt agaaaatgtt gaggttgag aatactcata ttatgattct 120
aagaatggag aaacttttga taagcaaatt ttatatcatt atccaatcct aaacgataag 180
ttaaaaatag gtaaattttg ctcaatagga ccagggtgaa ctattattat gaatggagca 240
aatcatagaa tggatggctc aacatatcca ttttaattat ttgtaatgg atgggagaaa 300
catatgccaa aattagatca actacctatt aagggggata caataatagg taatgatgta 360
tggataggaa aagatgttgt aattatgcca ggagtaaaaa tcggggatgg tgcaatagta 420
gctgctaatt ctgttttgtt aaaagatata gcgccataca tgttagctgg aggaaatcct 480
gctaacgaaa taaaacaaag atttgatcaa gatacaataa atcagctgct tgatataaaa 540
tggtggaatt ggccaataga cattattaat gagaatatag ataaaattct tgataatagc 600
atcattagag aagtcatatg gaaaaaatga 630

```

```

<210> SEQ ID NO 174
<211> LENGTH: 1440
<212> TYPE: DNA
<213> ORGANISM: Enterococcus faecalis

```

```
<400> SEQUENCE: 174
```

```

atgaatatag ttgaaaaatga aatatgtata agaactttaa tagatgatga ttttcctttg 60
atgttaaaat ggtaaactga tgaaagatga ttagaatttt atggtggtag agataaaaaa 120
tatacattag aatcattaaa aaaacattat acagagcctt ggaagatga agtttttaga 180
gtaattattg aatataacaa tgttcctatt ggatattgac aaatatataa aatgtatgat 240
gagttatata ctgattatca ttatccaaaa actgatgaga tagtctatgg tatggatcaa 300
tttataggag agccaaatta ttggagtaaa ggaattggta caagatatat taaattgatt 360
tttgaatttt tgaaaaaaga aagaaatgct aatgcagtta ttttagacc tcataaaaat 420
aatccaagag caataagggc ataccaaaaa tctggtttta gaattattga agatttgcca 480
gaacatgaat tacacgaggg caaaaaagaa gattgttatt taatggaata tagatatgat 540
gataatgcca caaatgttaa ggcaatgaaa tatttaattg agcattactt tgataatttc 600
aaagtagata gtattgaaat aatcggtagt ggttatgata gtgtggcata tttagttaat 660
aatgaataca tttttaaaac aaaatttagt actaataaga aaaaaggtta tgcaaaagaa 720
aaagcaatat ataatttttt aaatacaaat ttagaaacta atgtaaaaat tcctaattt 780
gaatattcgt atattagtga tgaattatct atactaggtt ataagaaat taaaggact 840
tttttaacac cagaaattta ttctactatg tcagaagaag aacaaaattt gttaaaacga 900

```

-continued

```

gatattgcca gttttttaag acaaatgcac ggttttagatt atacagatat tagtgaatgt   960
actattgata ataaacaaaa tgtattagaa gagtatatat tgttgcgtga aactatttat   1020
aatgatttaa ctgatataga aaaagattat atagaaagtt ttatggaaag actaaaatgca   1080
acaacagttt ttgagggtaa aaagtgttta tgccataatg attttagttg taatcatcta   1140
ttgttagatg gcaataatag attaactgga ataattgatt ttggagattc tggaattata   1200
gatgaatatt gtgattttat atacttactt gaagatagtg aagaagaaat aggaacaaat   1260
tttgagaag atatatataag aatgatgga aatatagata ttgagaaagc aaaagaatat   1320
caagatatag ttgaagaata ttatcctatt gaaactattg tttatggaat taaaaatatt   1380
aaacaggaat ttatcgaaaa tggtagaaaa gaaatttata aaaggactta taaagattga   1440

```

<210> SEQ ID NO 175

<211> LENGTH: 660

<212> TYPE: DNA

<213> ORGANISM: Staphylococcus aureus

<400> SEQUENCE: 175

```

ttgaatttaa acaatgacca tggacctgat cccgaaaata ttttaccgat aaaagggat   60
cggaatcttc aatttataaa acctactata acgaacgaaa acattttggt gggggaatat   120
tcttattatg atagtaagcg aggagaatcc tttgaagatc aagtcttata tcattatgaa   180
gtgattggag ataagttgat tataggaaga ttttgttcaa ttggtcccg aacaacattt   240
attatgaatg gtgcaaacca tcggatggat ggatcaacat atccttttca tctattcagg   300
atgggttggg agaagtatat gccttcttta aaagatcttc ccttgaaagg ggacattgaa   360
attgaaatg atgtatggat aggtagagat gtaaccatta tgcctggggg gaaaattggg   420
gacggggcaa tcattgctgc agaagctgtt gtcacaaaga atgttgctcc ctattctatt   480
gtcggtgtaa atcccttaaa atttataaga aaaaggtttt ctgatggagt tatcgaagaa   540
tggttagctt tacaatggtg gaatttagat atgaaaatta ttaatgaaa tcttcccttc   600
ataataaatg gagatatcga aatgtgtaag agaaaaagaa aacttctaga tgacacttga   660

```

<210> SEQ ID NO 176

<211> LENGTH: 1569

<212> TYPE: DNA

<213> ORGANISM: Staphylococcus aureus

<400> SEQUENCE: 176

```

atgaaaataa tgtagaggg acttaatata aaacattatg ttcaagatcg tttattgttg   60
aacataaatc gcctaaagat ttatcagaat gatcgtattg gtttaattgg taaaaatgga   120
agtggaaaaa caacgttact tcacatatta tataaaaaaa ttgtgcctga agaaggtatt   180
gtaaaacaat tttcacattg tgaacttatt cctcaattga agctcataga atcaactaaa   240
agtgtgggtg aagtaacacg aaactatatt cggcaagcgc ttgataaaaa tccagaactg   300
ctattagcag atgaaccaac aactaactta gataataact atatagaaaa attagaacag   360
gatttaaaaa attggcatg agcatttatt atagtttcac atgatcgcgc ttttttagat   420
aacttgtgta ctactatag ggaattgac gagggagaa taactgaata taaggggaat   480
tatagtaact atgttgaa caaagaatta gaaagacatc gagaagaatt agaatatgaa   540
aaatatgaaa aagaaaagaa acgattggaa aaagctataa atataaaaga acagaaagct   600

```

-continued

caacgagcaa ctaaaaaacc gaaaaactta agtttatctg aaggcaaaat aaaaggagca	660
aagccatact ttgcaggtaa gcaaagaag ttacgaaaaa ctgtaaaatc tctagaaacc	720
agactagaaa aacttgaaag cgtcgaaaag agaaacgaac ttcctccact taaaatggat	780
ttagtgaact tagaaagtgt aaaaaataga actataatac gtggtgaaga tgtctcgggt	840
acaattgaag gacgggtatt gtggaagca aaaagtttta gtattcgcgg aggagacaag	900
atggcaatta tcggatctaa tggtagacga aagacaacgt ttattaaaaa aattgtgcat	960
gggaatcctg gtatttcatt atcgccatct gtcaaaatcg gttattttag ccaaaaaata	1020
gatacattag aattagataa gagcatttta gaaaatgttc aatcttcttc acaacaaaat	1080
gaaactctta tcgaaactat tctagctaga atgcattttt tttagagatga tgtttataaa	1140
ccaataagtg tcttaagtgg tggagagcga gttaaagtag cactaactaa agtattctta	1200
agtgaagtta atacgttggg actagatgaa ccaacaaact ttcttgatat ggaagctata	1260
gaggcgtttg aatctttgtt aaaggaatat aatggcagta taatctttgt atctcacgat	1320
cgtaaattha tcgaaaaagt agccactcga ataatgacaa ttgataataa agaataaaa	1380
atatttgatg gcacatatga acaatttaaa caagctgaaa agccaacaag gaatattaaa	1440
gaagataaaa aacttttact tgagacaaaa attacagaag tactcagtcg attgagtatt	1500
gaaccttcgg aagaattaga acaagagttt caaaacttaa taaatgaaaa aagaaatttg	1560
gataaataa	1569

<210> SEQ ID NO 177

<211> LENGTH: 1467

<212> TYPE: DNA

<213> ORGANISM: Staphylococcus epidermidis

<400> SEQUENCE: 177

atggaacaat atacaattaa atttaaccaa atcaatcata aattgacaga tttacgatca	60
cttaacatcg atcatcttta tgcttaccaa tttgaaaaaa tagcacttat tgggggtaat	120
ggtactggta aaaccacatt actaaatatg attgctcaa aaacaaaacc agaactctgga	180
acagttgaaa cgaatggcga aattcaatat tttgaacagc ttaacatgga tgtggaaaat	240
gattttaaca cgtagacgg tagtttaatg agtgaactcc atatacctat gcatacaacc	300
gacagtatga gtggtgggta aaaagcaaaa tataaattac gtaatgtcat atcaaattat	360
agtcgatat tactttttaga tgaacctaca aatcacttgg ataaaattgg taaagattat	420
ctgaataata ttttaaaata ttactatggt actttaatta tagtaagtca cgatagagca	480
cttagagacc aaattgctga cacaatttgg gatatacaag aagatggcac aataagagtg	540
tttaaaggta attacacaca gtatcaaaat caatatgaac aagaacagtt agaacaacaa	600
cgtaaatatg aacagtatat aagtgaaaaa caaagattgt cccaagccag taaagctaaa	660
cgaaatcaag cgcaacaaat ggcacaagca tcatcaaac aaaaaataa aagtatagca	720
ccagatcgtt taagtgcac aaagaaaaa ggcacggttg agaaggctgc tcaaaaacaa	780
gctaagcata ttgaaaaaag aatggaacat ttggaagaag ttgaaaaacc acaaagtat	840
catgaattca attttccaca aaataaaatt tatgatatcc ataataatta tccaatcatt	900
gcacaaaatc taacattggt taaaggaagt caaaaactgc taacacaagt acgattccaa	960
ataccataty gcaaaaaatat agcgctcgta ggtgcaaatg gtgtaggtaa gacaacttta	1020

-continued

cttgaagcta ttaccacca aatagagga attgattgtt ctctaaagt gcaaatggca	1080
tactatcgtc aacttgctta tgaagacatg cgtgacgttt cattattgca atatttaatg	1140
gatgaaacgg attcatcaga atcattcagt agagctatgt taaataactt gggtttaa	1200
gaagcacttg agcgttcttg taatgttttg agtgggtggg aaagaacgaa attatcgta	1260
gcagtattat ttcaacgaa agcgaatag ttaatttttg atgaaccaac taatttttta	1320
gatattaaaa cattagaagc attagaaatg tttatgaata aatcctctg aatcattttg	1380
ttcatcac atgatacaag gtttgtaaa catgtatcag ataaaaatg ggaattaaca	1440
ggacaatcta ttcattgat aacttaa	1467

What is claimed is

1. A method using probes (fragments and/or oligonucleotides) and/or amplification primers which are specific, ubiquitous and sensitive for determining the presence and/or amount of nucleic acids from bacterial species selected from the group consisting of *Escherichia coli*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa*, *Proteus mirabilis*, *Streptococcus pneumoniae*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Enterococcus faecalis*, *Staphylococcus saprophyticus*, *Streptococcus pyogenes*, *Haemophilus influenzae* and *Moraxella catarrhalis* in a any sample suspected of containing said bacterial nucleic acid, wherein said bacterial nucleic acid or variant or part thereof comprises a selected target region hybridizable with said probes or primers; said method comprising the steps of contacting said sample with said probes or primers and detecting the presence and/or amount of hybridized probes and/or amplified products as an indication of the presence and/or amount of said bacterial species.

2. A method as defined in claim 1 further using probes (fragments and/or oligonucleotides) and/or amplification primers which are universal and sensitive for determining the presence and/or amount of nucleic acids from any bacteria from any sample suspected of containing said bacterial nucleic acid, wherein said bacterial nucleic acid or variant or part thereof comprises a selected target region hybridizable with said probes or primers; said method comprising the steps of contacting said sample with said probes or primers and detecting the presence and/or amount of hybridized probes and/or amplified products as an indication of the presence and/or amount of said any bacteria.

3. A method as defined in claim 1 further using probes (fragments and/or oligonucleotides) and/or amplification primers which are specific, ubiquitous and sensitive for determining the presence and/or amount of nucleic acids from an antibiotic resistance gene selected from the group consisting of bla_{tem}, Bla_{rob}, Bla_{shv}, aadB, aacC1, aacC2, aacC3, aacA4, mecA, vanA, vanH, vanX, satA, aacA-aphD, vat, vga, msrA, sul and int in any sample suspected of containing said bacterial nucleic acid, wherein said bacterial nucleic acid or variant or part thereof comprises a selected target region hybridizable with said probes or primers; said method comprising the steps of contacting said sample with said probes or primers and detecting the presence and/or amount of hybridized probes and/or amplified products as an indication of the presence and/or amount of said antibiotic resistance gene.

4. The method of any one of claims 1, 2 and 3 which is performed directly on a sample obtained from human patients, animals, environment or food.

5. The method of any one of claims 1, 2 and 3 which is performed directly on a sample consisting of one or more bacterial colonies.

6. The method of any one of claims 1 to 5, wherein the bacterial nucleic acid is amplified by a method selected from the group consisting of:

- polymerase chain reaction (PCR),
- ligase chain reaction,
- nucleic acid sequence-based amplification,
- self-sustained sequence replication,
- strand displacement amplification,
- branched DNA signal amplification,
- nested PCR, and
- multiplex PCR.

7. The method of claim 6 wherein said bacterial nucleic acid is amplified by PCR.

8. The method of claim 7 wherein the PCR protocol is modified to determine within one hour the presence of said bacterial nucleic acids by performing for each amplification cycle an annealing step of only one second at 55° C. and a denaturation step of only one second at 95° C. without any elongation step.

9. A method for the detection, identification and/or quantification of *Escherichia coli* directly from a test sample or from bacterial colonies, which comprises the following steps:

- depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic

acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO:3, SEQ ID NO: 4, SEQ ID NO: 5, SEQ ID NO: 6, SEQ ID NO: 7, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Escherichia coli*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Escherichia coli* in said test sample.

10. A method as defined in claim 9, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-227 nucleotides in length which sequence is comprised in SEQ ID NO: 3 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-278 nucleotides in length which sequence is comprised in SEQ ID NO: 4 or a complementary sequence thereof,
- 3) an oligonucleotide of 12-1596 nucleotides in length which sequence is comprised in SEQ ID NO: 5 or a complementary sequence thereof,
- 4) an oligonucleotide of 12-2703 nucleotides in length which sequence is comprised in SEQ ID NO: 6 or a complementary sequence thereof,
- 5) an oligonucleotide of 12-1391 nucleotides in length which sequence is comprised in SEQ ID NO: 7 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Escherichia coli*.

11. The method of claim 10, wherein the probe for detecting nucleic acid sequences from *Escherichia coli* is selected from the group consisting of SEQ ID NO: 44, SEQ ID NO: 45, SEQ ID NO: 48, SEQ ID NO: 49, SEQ ID NO: 50, SEQ ID NO: 51, SEQ ID NO: 52, SEQ ID NO: 53, SEQ ID NO: 54 and a sequence complementary thereof.

12. A method for detecting the presence and/or amount of *Escherichia coli* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Escherichia coli* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 3, SEQ ID NO: 4, SEQ ID NO: 5, SEQ ID NO: 6 and SEQ ID NO: 7;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Escherichia coli* in said test sample.

13. The method of claim 12, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 42 and SEQ ID NO: 43,
- b) SEQ ID NO: 46 and SEQ ID NO: 47,
- c) SEQ ID NO: 55 and SEQ ID NO: 56, and
- d) SEQ ID NO: 131 and SEQ ID NO: 132.

14. A method for the detection, identification and/or quantification of *Klebsiella pneumoniae* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO:8, SEQ ID NO: 9, SEQ ID NO: 10, SEQ ID NO: 11, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Klebsiella pneumoniae*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Klebsiella pneumoniae* in said test sample.

15. A method as defined in claim 14, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-238 nucleotides in length which sequence is comprised in SEQ ID NO: 8 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-385 nucleotides in length which sequence is comprised in SEQ ID NO: 9 or a complementary sequence thereof,
- 3) an oligonucleotide of 12-462 nucleotides in length which sequence is comprised in SEQ ID NO: 10 or a complementary sequence thereof,

- 4) an oligonucleotide of 12-730 nucleotides in length which sequence is comprised in SEQ ID NO: 11 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Klebsiella pneumoniae*.

16. The method of claim 15, wherein the probe for detecting nucleic acid sequences from *Klebsiella pneumoniae* is selected from the group consisting of SEQ ID NO: 57, SEQ ID NO: 58, SEQ ID NO: 59, SEQ ID NO: 60, SEQ ID NO: 63, SEQ ID NO: 64, SEQ ID NO: 65, SEQ ID NO: 66, SEQ ID NO: 69 and a sequence complementary thereof.

17. A method for detecting the presence and/or amount of *Klebsiella pneumoniae* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Klebsiella pneumoniae* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 8, SEQ ID NO: 9, SEQ ID NO: 10, and SEQ ID NO: 11;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Klebsiella pneumoniae* in said test sample.

18. The method of claim 17, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 61 and SEQ ID NO: 62,
- b) SEQ ID NO: 67 and SEQ ID NO: 68,
- c) SEQ ID NO: 135 and SEQ ID NO: 136, and
- d) SEQ ID NO: 137 and SEQ ID NO: 138.

19. A method for the detection, identification and/or quantification of *Proteus mirabilis* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 12, SEQ ID NO: 13,

SEQ ID NO: 14, SEQ ID NO: 15, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Proteus mirabilis*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Proteus mirabilis* in said test sample.

20. A method as defined in claim 19, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-225 nucleotides in length which sequence is comprised in SEQ ID NO: 12 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-402 nucleotides in length which sequence is comprised in SEQ ID NO: 13 or a complementary sequence thereof,
- 3) an oligonucleotide of 12-157 nucleotides in length which sequence is comprised in SEQ ID NO: 14 or a complementary sequence thereof,
- 4) an oligonucleotide of 12-1348 nucleotides in length which sequence is comprised in SEQ ID NO: 15 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Proteus mirabilis*.

21. The method of claim 20, wherein the probe for detecting nucleic acid sequences from *Proteus mirabilis* is selected from the group consisting of SEQ ID NO: 70, SEQ ID NO: 71, SEQ ID NO: 72, SEQ ID NO: 73, SEQ ID NO: 76, SEQ ID NO: 77, SEQ ID NO: 80, SEQ ID NO: 81, SEQ ID NO: 82 and a sequence complementary thereof.

22. A method for detecting the presence and/or amount of *Proteus mirabilis* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Proteus mirabilis* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 12, SEQ ID NO: 13, SEQ ID NO: 14, and SEQ ID NO: 15;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Proteus mirabilis* in said test sample.

23. The method of claim 22, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 74 and SEQ ID NO: 75, and
- b) SEQ ID NO: 133 and SEQ ID NO: 134.

24. A method for the detection, identification and/or quantification of *Staphylococcus saprophyticus* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 21, SEQ ID NO: 22, SEQ ID NO: 23, SEQ ID NO: 24, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Staphylococcus saprophyticus*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Staphylococcus saprophyticus* in said test sample.

25. A method as defined in claim 24, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-172 nucleotides in length which sequence is comprised in SEQ ID NO: 21 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-155 nucleotides in length which sequence is comprised in SEQ ID NO: 22 or a complementary sequence thereof,
- 3) an oligonucleotide of 12-145 nucleotides in length which sequence is comprised in SEQ ID NO: 23 or a complementary sequence thereof,
- 4) an oligonucleotide of 12-265 nucleotides in length which sequence is comprised in SEQ ID NO: 24 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Staphylococcus saprophyticus*.

26. The method of claim 25, wherein the probe for detecting nucleic acid sequences from *Staphylococcus saprophyticus* is selected from the group consisting of SEQ

ID NO: 96, SEQ ID NO: 97, SEQ ID NO: 100, SEQ ID NO: 101, SEQ ID NO: 102, SEQ ID NO: 103, SEQ ID NO: 104 and a sequence complementary thereof.

27. A method for detecting the presence and/or amount of *Staphylococcus saprophyticus* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Staphylococcus saprophyticus* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 21, SEQ ID NO: 22, SEQ ID NO: 23, and SEQ ID NO: 24;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Staphylococcus saprophyticus* in said test sample.

28. The method of claim 27, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 98 and SEQ ID NO: 99, and
- b) SEQ ID NO: 139 and SEQ ID NO: 140.

29. A method for the detection, identification and/or quantification of *Moraxella catarrhalis* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 28, SEQ ID NO: 29, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Moraxella catarrhalis*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said label-

ling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Moraxella catarrhalis* in said test sample.

30. A method as defined in claim 29, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-526 nucleotides in length which sequence is comprised in SEQ ID NO: 28 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-466 nucleotides in length which sequence is comprised in SEQ ID NO: 29 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Moraxella catarrhalis*.

31. The method of claim 30, wherein the probe for detecting nucleic acid sequences from *Moraxella catarrhalis* is selected from the group consisting of SEQ ID NO: 108, SEQ ID NO: 109, SEQ ID NO: 110, SEQ ID NO: 111, SEQ ID NO: 114, SEQ ID NO: 115, SEQ ID NO: 116, SEQ ID NO: 117 and a sequence complementary thereof.

32. A method for detecting the presence and/or amount of *Moraxella catarrhalis* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Moraxella catarrhalis* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 28 and SEQ ID NO: 29;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Moraxella catarrhalis* in said test sample.

33. The method of claim 32, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 112 and SEQ ID NO: 113,
- b) SEQ ID NO: 118 and SEQ ID NO: 119, and
- c) SEQ ID NO: 160 and SEQ ID NO: 119.

34. A method for the detection, identification and/or quantification of *Pseudomonas aeruginosa* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotide sequence is selected from the group consisting of SEQ ID NO: 16, SEQ ID NO: 17, SEQ ID NO: 18, SEQ ID NO: 19, SEQ ID NO: 20, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Pseudomonas aeruginosa*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Pseudomonas aeruginosa* in said test sample.

35. A method as defined in claim 34, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-2167 nucleotides in length which sequence is comprised in SEQ ID NO: 16 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-1872 nucleotides in length which sequence is comprised in SEQ ID NO: 17 or a complementary sequence thereof,
- 3) an oligonucleotide of 12-3451 nucleotides in length which sequence is comprised in SEQ ID NO: 18 or a complementary sequence thereof,
- 4) an oligonucleotide of 12-744 nucleotides in length which sequence is comprised in SEQ ID NO: 19 or a complementary sequence thereof,
- 5) an oligonucleotide of 12-2760 nucleotides in length which sequence is comprised in SEQ ID NO: 20 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Pseudomonas aeruginosa*.

36. The method of claim 35, wherein the probe for detecting nucleic acid sequences from *Pseudomonas aeruginosa* is selected from the group consisting of SEQ ID NO: 87, SEQ ID NO: 88, SEQ ID NO: 89, SEQ ID NO: 90, SEQ ID NO: 91, SEQ ID NO: 92, SEQ ID NO: 93, SEQ ID NO: 94, SEQ ID NO: 95 and a sequence complementary thereof.

37. A method for detecting the presence and/or amount of *Pseudomonas aeruginosa* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers

being capable of hybridizing selectively with one of the two complementary strands of *Pseudomonas aeruginosa* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 16, SEQ ID NO: 17, SEQ ID NO: 18, SEQ ID NO: 19 and SEQ ID NO: 20;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Pseudomonas aeruginosa* in said test sample.

38. The method of claim 37, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 83 and SEQ ID NO: 84, and
- b) SEQ ID NO: 85 and SEQ ID NO: 86.

39. A method for the detection, identification and/or quantification of *Staphylococcus epidermidis* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotide sequence is selected from the group consisting of SEQ ID NO: 36, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Staphylococcus epidermidis*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Staphylococcus epidermidis* in said test sample.

40. A method as defined in claim 39, wherein said probe is selected from the group consisting of an oligonucleotide of 12-705 nucleotides in length which sequence is comprised in SEQ ID NO: 36 and variants thereof which

specifically and ubiquitously anneal with strains and representatives of *Staphylococcus epidermidis*.

41. A method for detecting the presence and/or amount of *Staphylococcus epidermidis* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Staphylococcus epidermidis* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the following sequence: SEQ ID NO: 36;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Staphylococcus epidermidis* in said test sample.

42. The method of claim 41, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 145 and SEQ ID NO: 146, and
- b) SEQ ID NO: 147 and SEQ ID NO: 148.

43. A method for the detection, identification and/or quantification of *Staphylococcus aureus* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotide sequence is selected from the group consisting of SEQ ID NO: 37, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Staphylococcus aureus*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Staphylococcus aureus* in said test sample.

44. A method as defined in claim 43, wherein said probe is selected from the group consisting of an oligonucleotide of 12-442 nucleotides in length which sequence is comprised in SEQ ID NO: 37 and variants thereof which specifically and ubiquitously anneal with strains and representatives of *Staphylococcus aureus*.

45. A method for detecting the presence and/or amount of *Staphylococcus aureus* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Staphylococcus aureus* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the following sequence: SEQ ID NO: 37;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Staphylococcus aureus* in said test sample.

46. The method of claim 45, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 149 and SEQ ID NO: 150,
- b) SEQ ID NO: 149 and SEQ ID NO: 151, and
- c) SEQ ID NO: 152 and SEQ ID NO: 153.

47. A method for the detection, identification and/or quantification of *Haemophilus influenzae* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 25, SEQ ID NO: 26, SEQ ID NO: 27, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Haemophilus influenzae*, under conditions such that

the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Haemophilus influenzae* in said test sample.

48. A method as defined in claim 47, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-845 nucleotides in length which sequence is comprised in SEQ ID NO: 25 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-1598 nucleotides in length which sequence is comprised in SEQ ID NO: 26 or a complementary sequence thereof,
- 3) an oligonucleotide of 12-9100 nucleotides in length which sequence is comprised in SEQ ID NO: 27 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Haemophilus influenzae*.

49. The method of claim 48, wherein the probe for detecting nucleic acid sequences from *Haemophilus influenzae* is selected from the group consisting of SEQ ID NO: 105, SEQ ID NO: 106, SEQ ID NO: 107 and a sequence complementary thereof.

50. A method for detecting the presence and/or amount of *Haemophilus influenzae* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Haemophilus influenzae* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 25, SEQ ID NO: 26 and SEQ ID NO: 27;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Haemophilus influenzae* in said test sample.

51. The method of claim 50, wherein said at least one pair of primers comprises the following pair: SEQ ID NO: 154 and SEQ ID NO: 155.

52. A method for the detection, identification and/or quantification of *Streptococcus pneumoniae* directly from a test sample or from bacterial colonies, which comprises the following steps:

a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 30, SEQ ID NO: 31, SEQ ID NO: 34, SEQ ID NO: 35, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Streptococcus pneumoniae*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Streptococcus pneumoniae* in said test sample.

53. A method as defined in claim 52, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-631 nucleotides in length which sequence is comprised in SEQ ID NO: 30 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-3754 nucleotides in length which sequence is comprised in SEQ ID NO: 31 or a complementary sequence thereof,
- 3) an oligonucleotide of 12-841 nucleotides in length which sequence is comprised in SEQ ID NO: 34 or a complementary sequence thereof,
- 4) an oligonucleotide of 12-4500 nucleotides in length which sequence is comprised in SEQ ID NO: 35 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Streptococcus pneumoniae*.

54. The method of claim 53, wherein the probe for detecting nucleic acid sequences from *Streptococcus pneumoniae* is selected from the group consisting of SEQ ID NO: 120, SEQ ID NO: 121 and a sequence complementary thereof.

55. A method for detecting the presence and/or amount of *Streptococcus pneumoniae* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers

being capable of hybridizing selectively with one of the two complementary strands of *Streptococcus pneumoniae* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 30, SEQ ID NO: 31, SEQ ID NO: 34 and SEQ ID NO: 35;

b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Streptococcus pneumoniae* in said test sample.

56. The method of claim 55, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 78 and SEQ ID NO: 79,
- b) SEQ ID NO: 156 and SEQ ID NO: 157, and
- c) SEQ ID NO: 158 and SEQ ID NO: 159.

57. A method for the detection, identification and/or quantification of *Streptococcus pyogenes* directly from a test sample or from bacterial colonies, which comprises the following steps:

a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 32, SEQ ID NO: 33, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously anneals with strains or representatives of *Streptococcus pyogenes*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Streptococcus pyogenes* in said test sample.

58. A method as defined in claim 57, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-1337 nucleotides in length which sequence is comprised in SEQ ID NO: 32 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-1837 nucleotides in length which sequence is comprised in SEQ ID NO: 33 or a complementary sequence thereof, and

variants thereof which specifically and ubiquitously anneal with strains and representatives of *Streptococcus pyogenes*.

59. A method for detecting the presence and/or amount of *Streptococcus pyogenes* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Streptococcus pyogenes* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 32 and SEQ ID NO: 33;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Streptococcus pyogenes* in said test sample.

60. The method of claim 59, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 141 and SEQ ID NO: 142, and
- b) SEQ ID NO: 143 and SEQ ID NO: 144.

61. A method for the detection, identification and/or quantification of *Enterococcus faecalis* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 1, SEQ ID NO: 2, a sequence complementary thereof, a part thereof and a variant thereof, which specifically and ubiquitously

anneals with strains or representatives of *Enterococcus faecalis*, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of *Enterococcus faecalis* in said test sample.

62. A method as defined in claim 61, wherein said probe is selected from the group consisting of:

- 1) an oligonucleotide of 12-1817 nucleotides in length which sequence is comprised in SEQ ID NO: 1 or a complementary sequence thereof,
- 2) an oligonucleotide of 12-2275 nucleotides in length which sequence is comprised in SEQ ID NO: 2, and variants thereof which specifically and ubiquitously anneal with strains and representatives of *Enterococcus faecalis*.

63. A method for detecting the presence and/or amount of *Enterococcus faecalis* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of *Enterococcus faecalis* DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within one of the following sequences: SEQ ID NO: 1 and SEQ ID NO: 2;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of *Enterococcus faecalis* in said test sample.

64. The method of claim 63, wherein said at least one pair of primers is selected from the group consisting of:

- a) SEQ ID NO: 38 and SEQ ID NO: 39, and
- b) SEQ ID NO: 40 and SEQ ID NO: 41.

65. A method for the detection of the presence and/or amount of any bacterial species directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this

sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a universal probe which sequence is selected from the group consisting of SEQ ID NO: 122, SEQ ID NO: 123, SEQ ID NO: 124, SEQ ID NO: 125, SEQ ID NO: 128, SEQ ID NO: 129, SEQ ID NO: 130 and a sequence complementary thereof, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of the presence and/or amount of said any bacterial species in said test sample.

66. A method for detecting the presence and/or amount of any bacterial species in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing a pair of universal primers which sequence is defined in SEQ ID NO: 126 and SEQ ID NO: 127, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said any bacterial species DNA that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of the presence and/or amount of said any bacterial species in said test sample.

67. A method for evaluating a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene bla_{tem} (TEM-1) directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotide sequence is selected from the group consisting of SEQ ID NO: 161, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for a β -lactamase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene TEM-1.

68. A method as defined in claim 67, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 161.

69. A method for evaluating a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene bla_{tem} (TEM-1) in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for a β -lactamase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 161;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene TEM-1.

70. A method for evaluating a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene bla_{rob} (ROB-1) directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 162, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for a β -lactamase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene ROB-1.

71. A method as defined in claim 70, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 162.

72. A method for evaluating a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene bla_{rob} (ROB-1) in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for a β -lactamase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 162;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene ROB-1.

73. A method for evaluating a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene bla_{shv} (SHV-1) directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or
inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 163, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for a β -lactamase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene SHV-1.

74. A method as defined in claim 73, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 163.

75. A method for evaluating a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene bla_{shv} (SHV-1) in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for a β -lactamase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 163;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene SHV-1.

76. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene $aadB$ directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or
inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 164, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an aminoglycoside adenylyltransferase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aadB*.

77. A method as defined in claim 76, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 164.

78. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aadB* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an aminoglycoside adenylyltransferase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 164;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aadB*.

79. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC1* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or
inoculating said sample or said substantially homogenous population of bacteria isolated from this

sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 165, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC1*.

80. A method as defined in claim 79, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 165.

81. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC1* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 165;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC1*.

82. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC2* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 166, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC2*.

83. A method as defined in claim 82, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 166.

84. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC2* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 166;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC2*.

85. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC3* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 167, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC3*.

86. A method as defined in claim 85, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 167.

87. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC3* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 167;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacC3*.

88. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial anti-

biotic resistance gene *aacA4* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 168, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacA4*.

89. A method as defined in claim 88, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 168.

90. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacA4* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 168;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance

to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacA4*.

91. A method for evaluating a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene *mecA* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 169, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for a penicillin-binding protein, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene *mecA*.

92. A method as defined in claim 91, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 169.

93. A method for evaluating a bacterial resistance to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene *mecA* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for a penicillin-binding protein that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 169;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance

to β -lactam antibiotics mediated by the bacterial antibiotic resistance gene *mecA*.

94. A method for evaluating a bacterial resistance to vancomycin mediated by the bacterial antibiotic resistance genes *vanH*, *vanA* and *vanX* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 170, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance genes coding for vancomycin-resistance proteins, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to vancomycin mediated by the bacterial antibiotic resistance genes *vanH*, *vanA* and *vanX*.

95. A method as defined in claim 94, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 170.

96. A method for evaluating a bacterial resistance to vancomycin mediated by the bacterial antibiotic resistance genes *vanH*, *vanA* and *vanX* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance genes coding for vancomycin-resistance proteins that contain a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 170;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to vancomycin mediated by the bacterial antibiotic resistance genes *vanH*, *vanA* and *vanX*.

97. A method for evaluating a bacterial resistance to streptogramin A mediated by the bacterial antibiotic resistance gene *satA* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 173, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for a streptogramin A acetyltransferase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to streptogramin A mediated by the bacterial antibiotic resistance gene *satA*.

98. A method as defined in claim 97, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 173.

99. A method for evaluating a bacterial resistance to streptogramin A mediated by the bacterial antibiotic resistance gene *satA* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for streptogramin A acetyltransferase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 173;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to streptogramin A mediated by the bacterial antibiotic resistance gene *satA*.

100. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacA-aphD* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 174, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase-phosphotransferase under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacA-aphD*.

101. A method as defined in claim 100, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 174.

102. A method for evaluating a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacA-aphD* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an aminoglycoside acetyltransferase-phosphotransferase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 174;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and

- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to aminoglycoside antibiotics mediated by the bacterial antibiotic resistance gene *aacA-aphD*.

103. A method for evaluating a bacterial resistance to virginiamycin mediated by the bacterial antibiotic resistance gene *vat* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 175, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for a virginiamycin acetyltransferase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to virginiamycin mediated by the bacterial antibiotic resistance gene *vat*.

104. A method as defined in claim 103, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 175.

105. A method for evaluating a bacterial resistance to virginiamycin mediated by the bacterial antibiotic resistance gene *vat* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for a virginiamycin acetyltransferase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said

at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 175;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to virginiamycin mediated by the bacterial antibiotic resistance gene vat.

106. A method for evaluating a bacterial resistance to virginiamycin mediated by the bacterial antibiotic resistance gene vga directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 176, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an ATP-binding protein, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to virginiamycin mediated by the bacterial antibiotic resistance gene vga.

107. A method as defined in claim 106, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 176.

108. A method for evaluating a bacterial resistance to virginiamycin mediated by the bacterial antibiotic resistance gene vga in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an ATP-binding protein that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which

contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 176;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to virginiamycin mediated by the bacterial antibiotic resistance gene vga.

109. A method for evaluating a bacterial resistance to erythromycin mediated by the bacterial antibiotic resistance gene msrA directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 177, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an erythromycin resistance protein under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of a bacterial resistance to erythromycin mediated by the bacterial antibiotic resistance gene msrA.

110. A method as defined in claim 109, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 177.

111. A method for evaluating a bacterial resistance to erythromycin mediated by the bacterial antibiotic resistance gene msrA in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an erythromycin resistance protein that contains a target sequence, and the other of

said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 177;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of a bacterial resistance to erythromycin mediated by the bacterial antibiotic resistance gene *msrA*.

112. A method for evaluating potential bacterial resistance to β -lactams, aminoglycosides, chloramphenicol and/or trimethoprim mediated by the bacterial antibiotic resistance gene *int* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or

inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 171, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for an integrase, under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and
- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of potential bacterial resistance to β -lactams, aminoglycosides, chloramphenicol and/or trimethoprim mediated by the bacterial antibiotic resistance gene *int*.

113. A method as defined in claim 112, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 171.

114. A method for evaluating potential bacterial resistance to β -lactams, aminoglycosides, chloramphenicol and/or trimethoprim mediated by the bacterial antibiotic resistance gene *int* in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers

being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for an integrase that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 171;

- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of potential bacterial resistance to β -lactams, aminoglycosides, chloramphenicol and/or trimethoprim mediated by the bacterial antibiotic resistance gene *int*.

115. A method for evaluating potential bacterial resistance to β -lactams, aminoglycosides, chloramphenicol and/or trimethoprim mediated by the bacterial antibiotic resistance gene *int* directly from a test sample or from bacterial colonies, which comprises the following steps:

- a) depositing and fixing on an inert support or leaving in solution the bacterial DNA of the sample or of a substantially homogenous population of bacteria isolated from this sample, or
- inoculating said sample or said substantially homogenous population of bacteria isolated from this sample on an inert support, and lysing in situ said inoculated sample or isolated bacteria to release the bacterial DNA,

said bacterial DNA being in a substantially single stranded form;

- b) contacting said single stranded DNA with a probe, said probe comprising at least one single stranded nucleic acid which nucleotidic sequence is selected from the group consisting of SEQ ID NO: 172, a sequence complementary thereof, a part thereof and a variant thereof, which specifically anneals with said bacterial antibiotic resistance gene coding for a sulfonamide resistance protein under conditions such that the nucleic acid of said probe can selectively hybridize with said bacterial DNA, whereby a hybridization complex is formed, said complex being detected by labelling means, the label being present on said probe or the label being present on a first reactive member of said labelling means, said first reactive member reacting with a second reactive member present on said probe; and

- c) detecting the presence or the intensity of said label on said inert support or in said solution as an indication of potential bacterial resistance to β -lactams, aminoglycosides, chloramphenicol and/or trimethoprim mediated by the bacterial antibiotic resistance gene *int*.

116. A method as defined in claim 115, wherein said probe comprises an oligonucleotide of at least 12 nucleotides in length which hybridizes to SEQ ID NO: 172.

117. A method for evaluating potential bacterial resistance to β -lactams, aminoglycosides, chloramphenicol and/or tri-

methoprim mediated by the bacterial antibiotic resistance gene sul in a test sample which comprises the following steps:

- a) treating said sample with an aqueous solution containing at least one pair of oligonucleotide primers having at least 12 nucleotides in length, one of said primers being capable of hybridizing selectively with one of the two complementary strands of said bacterial antibiotic resistance gene coding for a sulfonamide resistance protein that contains a target sequence, and the other of said primers being capable of hybridizing with the other of said strands so as to form an extension product which contains the target sequence as a template, said at least one pair of primers being chosen from within the sequence defined in SEQ ID NO: 172;
- b) synthesizing an extension product of each of said primers which extension products contain the target sequence, and amplifying said target sequence, if any, to a detectable level; and
- c) detecting the presence and/or amount of said amplified target sequence as an indication of potential bacterial resistance to β -lactams, aminoglycosides, chloramphenicol and/or trimethoprim mediated by the bacterial antibiotic resistance gene sul.

118. A nucleic acid having the nucleotide sequence of any one of SEQ ID NOs: 1 to 37, SEQ ID NOs: 161 to 177, a part thereof and variants thereof which, when in single stranded form, ubiquitously and specifically hybridize with a target bacterial DNA as a probe or as a primer.

119. An oligonucleotide having a nucleotidic sequence of any one of SEQ ID NOs: 38 to 160.

120. A recombinant plasmid comprising a nucleic acid as defined in claim 118.

121. A recombinant host which has been transformed by a recombinant plasmid according to claim 120.

122. A recombinant host according to claim 121 wherein said host is *Escherichia coli*.

123. A diagnostic kit for the detection and/or quantification of the nucleic acids of any combination of the bacterial species defined in any one of claims **9, 14, 19, 24, 29, 34, 39, 43, 47, 52, 57** and **61**, comprising any combination of probes defined therein.

124. A diagnostic kit for the detection and/or quantification of the nucleic acids of any combination of the bacterial species defined in any one of claims **10, 11, 15, 16, 20, 21,**

25, 26, 30, 31, 35, 36, 40, 44, 48, 49, 53, 54, 58, 62 and **65**, comprising any combination of oligonucleotide probes defined therein.

125. A diagnostic kit for the detection and/or quantification of the nucleic acids of any combination of the bacterial species defined in any one of claims **12, 13, 17, 18, 22, 23, 27, 28, 32, 33, 37, 38, 41, 42, 45, 46, 50, 51, 55, 56, 59, 60, 63, 64** and **66** comprising any combination of primers defined therein.

126. A diagnostic kit for the detection and/or quantification of the nucleic acids of any combination of the bacterial resistance genes defined in any one of claims **67, 70, 73, 76, 79, 82, 85, 88, 91, 94, 97, 100, 103, 106** and **109** comprising any combination of probes defined therein.

127. A diagnostic kit for the detection and/or quantification of the nucleic acids of any combination of the bacterial resistance genes defined in any one of claims **68, 71, 74, 77, 80, 83, 86, 89, 92, 95, 98, 101, 104, 107** and **110** comprising any combination of oligonucleotide probes defined therein.

128. A diagnostic kit for the detection and/or quantification of the nucleic acids of any combination of the bacterial resistance genes defined in any one of claims **69, 72, 75, 78, 81, 84, 87, 90, 93, 96, 99, 102, 105, 108** and **111** comprising any combination of primers defined therein.

129. A diagnostic kit for the simultaneous detection and quantification of nucleic acids of any combination of the bacterial species defined in claim 123, comprising any combination of the bacterial probes defined therein and any combination of the probes to the antibiotic resistance genes defined in any one of SEQ ID NOs: 161 to 177 in whole or in part.

130. A diagnostic kit for the simultaneous detection and quantification of nucleic acids of any combination of the bacterial species defined in claim 124, comprising any combination of the bacterial oligonucleotide probes defined therein and any combination of oligonucleotide probes that hybridize to the antibiotic resistance genes defined in any one of SEQ ID NOs: 161 to 177.

131. A diagnostic kit for the simultaneous detection and quantification of nucleic acids of any combination of the bacterial species defined in claim 125, comprising any combination of the primers defined therein and any combination of primers that anneal to the antibiotic resistance genes defined in any one of SEQ ID NOs: 161 to 177.

* * * * *