



Office de la Propriété
Intellectuelle
du Canada

Un organisme
d'Industrie Canada

Canadian
Intellectual Property
Office

An agency of
Industry Canada

CA 2464874 C 2013/10/22

(11)(21) **2 464 874**

(12) **BREVET CANADIEN
CANADIAN PATENT**

(13) **C**

(22) Date de dépôt/Filing Date: 2004/04/30

(41) Mise à la disp. pub./Open to Public Insp.: 2004/11/01

(45) Date de délivrance/Issue Date: 2013/10/22

(30) Priorité/Priority: 2003/05/01 (US10/427,243)

(51) Cl.Int./Int.Cl. *C12Q 1/68* (2006.01),
A61B 10/00 (2006.01), *C12Q 1/04* (2006.01),
G01N 33/574 (2006.01)

(72) Inventeurs/Inventors:
ATKINS, DAVID, US;
BACKUS, JOHN, US;
BELLY, ROBERT, US;
ROSEN, STEVEN, US;
WHITE, ROBERT, US

(73) Propriétaire/Owner:
VERIDEX, LLC, US

(74) Agent: NORTON ROSE FULBRIGHT CANADA
LLP/S.E.N.C.R.L., S.R.L.

(54) Titre : ESSAI PEROPERATOIRE DU GANGLION LYMPHATIQUE

(54) Title: INTRAOPERATIVE LYMPH NODE ASSAY

(57) **Abrégé/Abstract:**

In an intraoperative lymph node assay for micrometastases, a surgeon identifies a sentinel lymph node during surgery according to known methods. Sentinel nodes are removed and prepared for the extraction of nucleic acid. Nucleic acid (e.g., RNA) is then rapidly extracted from the sentinel nodes. Markers indicative of micrometastases, if present, are then amplified and detected.



Docket No: CDS-5009

Abstract of the Disclosure

INTRAOPERATIVE LYMPH NODE ASSAY

In an intraoperative lymph node assay for micrometastases, a surgeon identifies a
5 sentinel lymph node during surgery according to known methods. Sentinel nodes
are removed and prepared for the extraction of nucleic acid. Nucleic acid (e.g.,
RNA) is then rapidly extracted from the sentinel nodes. Markers indicative of
micrometastases, if present, are then amplified and detected.

Docket No: CDS-5009

INTRAOPERATIVE LYMPH NODE ASSAY

BACKGROUND

5 The invention relates to the field of molecular diagnostics.

Lymph node involvement is the strongest prognostic factor in many solid tumors, and detection of lymph node micrometastases is of great interest to pathologists and surgeons. Current lymph node evaluation involves microscopic examination of H&E-stained tissue sections and suffers from three major limitations: (a) single
10 tumor cells, or small foci of cells, are easily missed; (b) the result is not rapidly available, meaning that any positive result in a sentinel lymph node procedure requires a second surgery for removal of axillary lymph nodes and (c) only one or two tissue sections are studied, and thus the vast majority of each node is left unexamined. Serial sectioning can help overcome the issue of sampling error, and
15 immunohistochemistry (IHC) can help identify individual tumor cells. The combination of these methods, however, is too costly and time consuming for routine analysis and is limited to special cases such as sentinel lymph node examination.

Surgical decisions are often based on intra-operative frozen section analysis of
20 lymph nodes; however, the sensitivity of these methods is relatively poor, ranging from 50-70% relative to standard H&E pathology, leading to an unacceptably high rate of second surgeries. The five-year survival of Stage 0 and I breast cancer patients who do not have lymph node involvement are 92% and 87%, respectively. On the other hand, the five-year survival of later stage breast cancer patients who do
25 have lymph node involvement decrease significantly. For example, the survival of Stage II breast cancer is only 75%, Stage III 46%, and Stage IV 13%. Although node negative breast cancer patients have improved survival, 20-30% of histologically

Docket No: CDS-5009

node negative patients suffer disease recurrence. This is most likely due to the limitations of current techniques in the detection of micrometastases including issues related to node sampling and poor sensitivity for detecting individual tumor cells or small tumor foci.

- 5 In addition to a need for more accurate and sensitive detection of metastases in breast cancer nodes, there is a similar need for a rapid and sensitive test in melanoma nodes as well as in the more accurate detection of surgical margins, particularly in the intraoperative setting.

- 10 The polymerase chain reaction (PCR) is a powerful tool in the field of molecular biology that could be useful in this setting. This technique allows for replicating/amplifying small amounts of nucleic acid fragments into quantities that can be analyzed in a meaningful way. Furthermore with the development of real-time quantitative RT-PCR (Q-RT-PCR), this technology has become more reliable as well as amenable to automation. Q-RT-PCR is less subject to contamination and
15 provides quantitation of gene expression. Such quantitation could be applied for the detection of micrometastases in intraoperative lymph node assays. PCR in molecular diagnostics, despite its advantages, has several limitations that make it difficult to apply in typical clinical diagnostic setting, particularly in the intraoperative setting.

- One such limitation is the time it typically takes to perform PCR diagnoses. Typical
20 PCR reactions take hours, not minutes. Decreasing the time it takes to carry out a PCR reaction is necessary if the technique is to be useful intraoperatively. Further, although Q-RT-PCR can provide quantitative results, to date there have been no known cutoff values for distinguishing positive from negative results based on such technology nor has it been clear which nucleic acid fragments are best detected and
25 correlated to the presence of a micrometastasis. Other methods for the amplification and detection of nucleic acid fragments exist as well and each suffers from similar problems.

Docket No: CDS-5009

An intraoperative molecular lymph node assay that overcomes the existing difficulties would be well accepted by the medical community.

SUMMARY OF THE INVENTION

The invention is an intraoperative lymph node assay for micrometastases. A
5 surgeon identifies a sentinel lymph node during surgery according to known methods. Sentinel nodes are removed and prepared as described below. Nucleic acid (e.g., RNA) is then rapidly extracted from the sentinel nodes. Markers indicative of micrometastases, if present, are then amplified and detected. The surgeon then takes action based upon the outcome of the detection of such Markers.

10 The assay is most preferably based on a rapid PCR and RT-PCR method that permits execution of a complete PCR reaction in minutes thereby permitting use of PCR in intraoperative diagnoses for detecting micrometastases in sentinel lymph nodes.

In one aspect of the invention, an intraoperative molecular diagnostic method
15 includes the following steps conducted during the course of a surgical procedure: obtaining a tissue sample from a patient; extracting RNA from the sample, analyzing the sample by nucleic acid amplification and detection; and determining if the presence of one or more gene expression Markers exceeds a cut-off value. The cut-off value can be an absolute value or a value relative to the expression of a control
20 gene.

In another aspect of the invention, the gene expression Markers are nucleic acid fragments specific for particular tissue such as breast tissue.

In yet another aspect of the invention, the gene expression Markers are nucleic acid fragments indicative of malignancy.

In yet another aspect of the invention lymph nodes samples are homogenized and then diluted to assure that sufficient RNA is extracted despite an expedited extraction method.

In yet another aspect of the invention the Markers are those of mammaglobin (Seq. ID No. 17) and either PIP (Seq. ID No. 18), B305D (particularly, isoform C, Seq. ID No. 14), B726 (Seq. ID No. 15), GABA. (Seq. ID No. 16) or O8E (Seq. ID No. __) in the case of breast cancer diagnostics and tyrosindse (Seq. ID No. 22) and MART 1 (Seq. ID No. 21) in the case of melanoma.

10 In yet another aspect of the invention, micrometastases are detected intraoperatively by a method that includes the steps of: obtaining RNA from a sentinel lymph node; performing a quantitative RT-PCR method specific to one or more genes of interest and determining if the presence of the Marker exceeds a predetermined cut-off. The cut-off value can be an absolute value or a value relative to the expression of a control
15 gene.

In a yet further embodiment of the invention, kits contain reagents for conducting intraoperative lymph node assay for micrometastases.

20 In yet another aspect, there is provided a method of conducting an intra-operative molecular diagnostic assay comprising the steps of: analyzing a lymph node tissue sample from a patient by extracting RNA from the lymph node; amplifying and detecting the RNA; and determining if the presence of more than one Marker exceeds a cut-off value, wherein said RNA is extracted in less than 8 minutes and wherein the method is
25 conducted within a period of no more than 35 minutes.

In yet another aspect, there is provided a method of conducting an intra-operative molecular diagnostic assay for use in detecting melanoma metastasis comprising the steps of: analyzing a lymph node tissue sample from a patient by extracting RNA from
30 the lymph node; amplifying and detecting the RNA; and determining if the presence of more than one Marker exceeds a cut-off value, wherein said RNA is extracted in less

than 8 minutes and wherein the method is conducted within a period of no more than 35 minutes.

5 **DETAILED DESCRIPTION**

Intraoperative methods for cancer diagnostics are presented. These methods employ a rapid method of extracting nucleic acids from a lymph node and a method of amplifying and detecting nucleic acid fragments indicative of metastasis (such fragments are
10 referred to herein as “Markers”).

An important aspect of the inventive methods is the rapid amplification and detection of Markers indicative of the expression of certain genes. Provided that such methods can be conducted within a period acceptable for an intraoperative assay (i.e., no more than about
15 35 minutes), any reliable, sensitive, and specific

method can be used. This includes PCR methods, Rolling Circle Amplification methods (RCA), Ligase Chain Reaction methods (LCR), Strand Displacement Amplification methods (SDA), Nucleic Acid Sequence Based Amplification methods (NASBA), and others. The rapid molecular diagnostics involved are most
5 preferably quantitative PCR methods, including QRT-PCR.

Irrespective of the amplification method employed, it is important to adequately sample the tissue used to conduct the assay. This includes proper excision and processing of the sentinel lymph node as well as extraction of RNA from it. Once
10 obtained, it is important to process the nodes properly so that any cancerous cells present are detected.

A variety of techniques are available for extracting nucleic acids from tissue samples. Standard practice in each case is time consuming and can be difficult even when using a commercially available kit designed for this purpose. Typical
15 commercially available nucleic acid extraction kits take at least 15 minutes to extract the nucleic acid. In the methods of the instant invention, nucleic acid is extracted in less than 8 minutes and preferably less than 6 minutes. These rapid extraction methods are the subject of a patent application by R. Belly, J. Toner, and J. Backus of equal date and entitled "Rapid Extraction Of RNA From Cells And Tissues".

The successful isolation of intact RNA generally involves four steps: effective
20 disruption of cells or tissue, denaturation of nucleoprotein complexes, inactivation of endogenous ribonuclease (RNAase) and removal of contaminating DNA and protein. The disruptive and protective properties of guanidinium thiocyanate (GTC) and B-mercaptoethanol to inactivate the ribonucleases present in cell extracts make them preferred reagents for the first step. When used in conjunction with a surfactant
25 such as sodium dodecylsulfate (SDS), disruption of nucleoprotein complexes is achieved allowing the RNA to be released into solution and isolated free of protein. Dilution of cell extracts in the presence of high concentrations of GTC causes

Docket No: CDS-5009

selective precipitation of cellular proteins to occur while RNA remains in solution. Centrifugation can clear the lysate of precipitated proteins and cellular DNA and is preferably performed through a column. Such columns also shear DNA and reduce the viscosity of the sample. RNA purification is preferably conducted on a spin
5 column containing silica or other material. Manual cell and tissue disruption can be by means of a disposable tissue grinder as described in US Patent 4,715,545. Homogenization time is within 1 to 2 minute and is more preferably 30-45sec. The sample can then processed with a shredding column (eg., QIAshredder, QIAGEN Inc., Valencia, CA, or suitable substitute) or with an RNA processing device such as
10 the PCR Tissue Homogenization Kit commercially available from Omni International (Warrenton, VA) to reduce its viscosity. RNA is precipitated out via the spin column as described above and centrifugation times are no greater than 30 sec. When using commercial RNA extraction kits such as those available from Qiagen, Inc., filtration is used instead of centrifugation for all steps except for the
15 column drying step. Typically, the sample is diluted with an equal volume of 70% ethanol prior to application on the column. After washes by filtration, the column is dried by centrifugation, and RNA is eluted in RNAase free water. The RNA is selectively precipitated out of solution with ethanol and bound to a substrate (preferably, a silica-containing membrane or filter). The binding of RNA to the
20 substrate occurs rapidly due to the disruption of the water molecules by the chaotropic salts, thus favoring absorption of nucleic acids to the silica. The bound total RNA is further purified from contaminating salts, proteins and cellular impurities by simple washing steps. Finally, the total RNA is eluted from the membrane by the addition of nuclease-free water. The total time of this rapid
25 protocol is less than 8 minutes and preferably less than 6 min.

In summary the rapid RNA extraction method involves the following steps:

- (a) obtaining a sample containing cells from the biological system,

Docket No: CDS-5009

- (b) optionally, removing from the sample, cells without RNA of interest to produce a working sample,
 - (c) lysing the cells containing RNA that is of interest and producing a homogenate of them,
 - 5 (d) optionally, diluting the homogenate,
 - (e) contacting the wetted, homogenized working sample with a substrate containing, or to which is affixed, a material to which RNA binds,
 - (f) allowing the sample to bind to the substrate,
 - (g) removing contaminants and interferents,
 - 10 (h) drying the substrate, and
 - (i) eluting RNA from the substrate;
- in instances in which centrifugation is used, it may occur after steps g, h, or i and vacuum/filtration is preferably applied in extraction steps.

15 The reagents involved in this rapid extraction process are:

- Lysis/Binding buffer (preferably, 4.5M guanidinium-HCl, 100mM sodium phosphate),
- Wash buffer I (preferably, 37% ethanol in 5M guanidine-HCL, 20mM Tris-HCL),
- 20 Wash buffer II (preferably, 80% ethanol in 20mM NaCl, 2mM Tris-HCl),
- Elution buffer, and
- Nuclease-free sterile double distilled water.

25 Since the distribution of cancer cells in nodes is non-uniform, it is preferable that multiple sections of the node be sampled. Optionally, one or more nodes may also be examined based on pathology. One method for accomplishing both a molecular based test and an examination of the same node sample by pathology is to section the node into at least four sections with one outer and inner section used for pathology, and one outer and inner section for used for molecular testing. As the

Docket No: CDS-5009

distribution of metastases and micrometastases in tissues is not uniform in nodes or other tissues, a sufficiently large sample should be obtained so that metastases will not be missed. One approach to this sampling issue in the present method is to homogenize a large tissue sample, and subsequently perform a dilution of the well-mixed homogenized sample to be used in subsequent molecular testing.

A typical PCR reaction includes multiple amplification steps, or cycles that selectively amplify target nucleic acid species. A typical PCR reaction includes three steps: a denaturing step in which a target nucleic acid is denatured; an annealing step in which a set of PCR primers (forward and backward primers) anneal to complementary DNA strands; and an elongation step in which a thermostable DNA polymerase elongates the primers. By repeating this step multiple times, a DNA fragment is amplified to produce an amplicon, corresponding to the target DNA sequence. Typical PCR reactions include 20 or more cycles of denaturation, annealing and elongation. In many cases, the annealing and elongation steps can be performed concurrently, in which case the cycle contains only two steps.

In the inventive method, employing RT-PCR, the RT-PCR amplification reaction is conducted in a time suitable for intraoperative diagnosis, the lengths of each of these steps can be in the seconds range, rather than minutes. Specifically, with certain new thermal cyclers being capable of generating a thermal ramp rate of at least about 5°C per second, RT-PCR amplifications in 30 minutes or less are used. More preferably, amplifications are conducted in less than 25 minutes. With this in mind, the following times provided for each step of the PCR cycle does not include ramp times. The denaturation step may be conducted for times of 10 seconds or less. In fact, some thermal cyclers have settings for "0 seconds" which may be the optimal duration of the denaturation step. That is, it is enough that the thermal cycler reaches the denaturation temperature. The annealing and elongation steps are most preferably less than 10 seconds each, and when conducted at the same temperature,

the combination annealing/elongation step may be less than 10 seconds. Some homogeneous probe detection methods, however, may require a separate step for elongation to maximize rapid assay performance. In order to minimize both the total amplification time and the formation of non-specific side reactions, annealing
5 temperatures are typically above 50 C. More preferably annealing temperatures are above 55 C.

A single combined reaction for RT-PCR, with no experimenter intervention, is desirable for several reasons: (1) decreased risk of experimenter error, (2) decreased risk of target or product contamination and (3) increased assay speed. The reaction
10 can consist of either one or two polymerases. In the case of two polymerases, one of these enzymes is typically an RNA-based DNA polymerase (reverse transcriptase) and one is a thermostable DNA-based DNA polymerase. To maximize assay performance, it is preferable to employ a form of "hot start" technology for both of these enzymatic functions. US Patents 5,411,876 and 5,985,619 provide examples
15 of different "hot start" approaches.

Preferred methods include the use of one or more thermoactivation methods that sequester one or more of the components required for efficient DNA polymerization. US Patents 5,550,044 and 5,413,924 describe methods for preparing reagents for use in such methods. US Patent 6,403,341 describes a sequestering approach that
20 involves chemical alteration of one of the PCR reagent components.

In the most preferred embodiment, both RNA- and DNA-dependent polymerase activities reside in a single enzyme. Other components that are required for efficient amplification include nucleoside triphosphates, divalent salts and buffer components. In some instance, non-specific
25 nucleic acid and enzyme stabilizers may be beneficial.

The specificity of any given amplification-based molecular diagnostic relies heavily, but not exclusively, on the identity of the primer sets. The primer sets are pairs of forward and reverse oligonucleotide primers that anneal to a target DNA sequence to

Docket No: CDS-5009

permit amplification of the target sequence, thereby producing a target sequence-specific amplicon. The primers must be capable of amplifying Markers of the disease state of interest. In the case of the instant invention, these Markers are directed to breast cancer and melanoma.

- 5 In the case of breast cancer, the inventive method involves the amplification of a tissue marker specific for either breast tissue or breast cancer tissue. A combination of at least two Markers are used such that clinically significant and reliable detection of breast/and or cancer cells in lymph nodes is detected when present. Preferably, the Markers are amplified and detected in a single reaction vessel at the same time
10 (i.e., they are multiplexed). Most preferably, the primer sets are complimentary to nucleic acid fragments specific to those Markers. The Markers include mammaglobin (Seq. ID No. 17) and one or more (preferably one) of the following: B305D (Seq. ID No. 14), PIP (Seq. ID No. 18), B726 (Seq ID No. 15), GABA (Seq ID No. 16) or O8E (Seq ID No. ____). The combination of a tissue specific marker
15 and a cancer specific marker provide sensitivity and specificity that exceeds 90 % and 95 % respectively. Some Markers exist in various isoforms with certain of the isoforms being more specific for one tissue or cancer than others. In the case of B305D, the most preferred isoform is B305D isoform C . It is also the most preferred Marker in combination with the mammaglobin Marker.
- 20 In the case of melanoma, the inventive method involves the amplification of markers specific to melanoma cancer. Thus, the primer sets are complimentary to nucleic acid fragments specific to those markers. The Markers include one or more (preferably two or more) of the following: Tyrosinase (Seq. ID No. 22) and MART 1 (Seq. ID No. 21). These Markers provide sensitivity and specificity that exceeds
25 90 %.

The reaction must also contain some means of detection of a specific signal. This is preferably accomplished through the use of a reagent that detects a region of DNA sequence derived from polymerization of the target sequence of interest. Preferred

reagents for detection give a measurable signal differential when bound to a specific nucleic acid sequence of interest. Often, these methods involve nucleic acid probes that give increased fluorescence when bound to the sequence of interest. The progress of the PCR reactions of the inventive method are typically monitored by
5 analyzing the relative rates of amplicon production for each PCR primer set. Monitoring amplicons production may be achieved by a number of detection reagents and methods, including without limitation, fluorescent primers, fluorogenic probes and fluorescent dyes that bind double-stranded DNA, molecular beacons, Scorpions, and others. A common method of monitoring a PCR reaction employs a
10 fluorescent 5'nuclease assay. This method exploits the 5'nuclease activity of certain thermostable DNA polymerases (such as Taq or Tfl DNA polymerases) to cleave an oligomeric probe during the PCR process. The oligomer is selected to anneal to the amplified target sequence under elongation conditions. The probe typically has a fluorescent reporter on its 5'end and a fluorescent quencher of the reporter at the 3'
15 end. So long as the oligomer is intact, the fluorescent signal from the reporter is quenched. However, when the oligomer is digested during the elongation process, the fluorescent reporter is no longer in proximity to the quencher. The relative accumulation of free fluorescent reporter for a given amplicon may be compared to the accumulation of the same amplicons for a control sample and/or to that of a
20 control gene, such as, without limitation, β -Actin and PBDG (porphobilinogen deaminase) to determine the relative abundance of a given cDNA product of a given RNA in a RNA population. Products and reagents for the fluorescent 5'nuclease assay are readily available commercially, for instance from Applied Biosystems.

The preferred detection reagents are commonly referred to as "Scorpions" and are
25 described in US Patents 6,326,145 and 5,525,494.

These reagents include one or more molecules comprising a tailed primer and an integrated signaling system. The primer has a template binding region and a tail comprising a linker and a target binding region. The target binding region in the tail hybridizes to complementary sequence in an extension product of the primer.

Docket No: CDS-5009

This target specific hybridization event is coupled to a signaling system wherein hybridization leads to a detectable change. In PCR reactions the target binding region and the tail region are advantageously arranged such that the tail region remains single stranded, ie. uncopied. Thus the tail region is non-amplifiable in the
5 PCR amplification products. The linker comprises a blocking moiety which prevents polymerase mediated chain extension on the primer template.

Equipment and software also are readily available for controlling and monitoring amplicon accumulation in PCR and QRT-PCR including the Smart Cycler thermocycler commercially available from Cepheid of Sunnyvale, California, and
10 the ABI Prism 7700 Sequence Detection System, commercially available from Applied Biosystems.

In the preferred RT-PCR reactions, the amounts of certain reverse transcriptase and the PCR reaction components are atypical in order to take advantage of the faster ramp times of some thermal cyclers. Specifically, the primer concentrations are very
15 high.

Typical gene-specific primer concentrations for reverse transcriptase reactions are less than about 20 nM. To achieve a rapid reverse transcriptase reaction on the order of one to two minutes, the reverse transcriptase primer concentration was raised to greater than 20 nM, preferably at least about 50 nM, and typically about 100 nM.
20 Standard PCR primer concentrations range from 100 nM to 300 nM. Higher concentrations may be used in standard PCR reactions to compensate for T_m variations. However, for purposes herein, the referenced primer concentrations are for circumstances where no T_m compensation is needed. Proportionately higher concentrations of primers may be empirically determined and used if T_m
25 compensation is necessary or desired. To achieve rapid PCR reactions, the PCR primer concentrations typically are greater than 250 nM, preferably greater than about 300 nM and typically about 500 nM.

Docket No: CDS-5009

Commercially used diagnostics also preferably employ the use of one or more internal positive controls that confirms the operation of a particular amplification reaction for a negative result. Potential causes of false negative results that must be controlled for in an RT-PCR reaction include: inadequate RNA quantity, degradation of RNA, inhibition of RT and/or PCR, and experimenter error. In the case of gene expression assays, it is preferable to utilize a gene that is constitutively expressed in the tissue of interest. PBGD (Seq. ID No. 20) is a gene that is commonly used as an internal control due to several factors: it contains no known pseudogenes in humans, it is constitutively expressed in human tissues, and it is expressed at a relatively low level and therefore is less likely to cause inhibition of the amplification of target sequences of interest. Use of PBGD as a control minimizes or eliminates reporting erroneous results arising from all potential sources of false negative results.

In the commercialization of the described methods for QRT-PCR certain kits for detection of specific nucleic acids are particularly useful. In one embodiment, the kit includes reagents for amplifying and detecting Markers. Optionally, the kit includes sample preparation reagents and or articles (e.g., tubes) to extract nucleic acids from lymph node tissue. The kits may also include articles to minimize the risk of sample contamination (e.g., disposable scalpel and surface for lymph node dissection and preparation).

In a preferred kit, reagents necessary for the one-tube QRT-PCR process described above are included such as reverse transcriptase, a reverse transcriptase primer, a corresponding PCR primer set (preferably for Markers and controls), a thermostable DNA polymerase, such as Taq polymerase, and a suitable detection reagent(s), such as, without limitation, a scorpion probe, a probe for a fluorescent 5'nuclease assay, a molecular beacon probe, a single dye primer or a fluorescent dye specific to double-stranded DNA, such as ethidium bromide. The primers are preferably in quantities that yield the high concentrations described above. Thermostable DNA polymerases

are commonly and commercially available from a variety of manufacturers. Additional materials in the kit may include: suitable reaction tubes or vials, a barrier composition, typically a wax bead, optionally including magnesium; reaction mixtures (typically 10X) for the reverse transcriptase and the PCR stages, including
 5 necessary buffers and reagents such as dNTPs; nuclease-or RNase-free water; RNase inhibitor; control nucleic acid (s) and/or any additional buffers, compounds, co-factors, ionic constituents, proteins and enzymes, polymers, and the like that may be used in reverse transcriptase and/or PCR stages of QRT-PCR reactions. Optionally, the kits include nucleic acid extraction reagents and materials.

10 The following non-limiting examples help to further describe the invention.

Examples

Real-time PCR

Examples in the present invention are based on the use of real-time PCR. In real-time PCR the products of the polymerase chain reaction are monitored in real-time
 15 during the exponential phase of PCR rather than by an end-point measurement. Hence, the quantification of DNA and RNA is much more precise and reproducible. Fluorescence values are recorded during every cycle and represent the amount of product amplified to that point in the amplification reaction. The more templates present at the beginning of the reaction, the fewer number of cycles it takes to reach
 20 a point in which the fluorescent signal is first recorded as statistically significant above background, which is the definition of the (Ct) values. The concept of the threshold cycle (Ct) allows for accurate and reproducible quantification using fluorescence based RT-PCR. Homogeneous detection of PCR products are preferably performed based on: (a) double- stranded DNA binding dyes (e.g., SYBR
 25 Green), (b) fluorogenic probes (e.g., Taq ManTM probes, Molecular Beacons), and (c) direct labeled primers (e.g., Amplifluor primers).

Example 1. A comparison of RNA extracted by the Rapid Method and by prior art Method based on real-time PCR with Taqman assays.

A total of 16 H&E positive nodes and 15 H&E negative breast node samples were purchased from Genomics Collaborative. Two 30 mg pieces of tissue were cut from
5 each node.

Part I (Prior Art). One 30mg piece of tissue was processed as follows: 600ul of Buffer RLT was added and the tissue sample was homogenized manually for 20-40 sec by means of a disposable tissue grinder (cat 15704-126, VWR Scientific, West
10 Chester, PA). The tube was centrifuged for 3min at maximum speed in EppendorfTM model 5415C microcentrifuge. The supernatant fluid was transferred to a new tube, and 1 volume of 70% ethanol was added.

Sample (700μL) was applied to an RNeasy mini column placed on the QIAvac 24
15 vacuum manifold and the vacuum was applied. Buffer RW1 (700 μL) and allowed to filter through the column. Buffer RPE (500 μL) was pipetted onto the column and allowed to filter through. Another 500 μL of Buffer RPE was added to the column, and allowed to filter through. The RNeasy column was placed in a 2 ml collection tube and centrifuged in a microfuge at maximum speed for 1min. The RNA was then
20 eluted from the column, by transferring the RNeasy colum to a new 1.5 ml collection tube, adding 50ul of RNAase-free water, and centrifuging for 1min at 8000.X g.

Part II (Rapid Extraction). The second 30 mg piece of breast node tissue was
25 processed by a rapid protocol as follows: (1) The tissue piece was added to 600 μL of Buffer RLT and homogenized manually for 20-40 sec by means of a disposable tissue grinder. (2) The homogenate was centrifuged through a QIAshredder column for 30sec at maximum speed. (3) 1 volume of 70% ethanol was added to the lysate and mixed by pipetting. (4) The sample was then applied to an RNeasy mini column

placed on the QIAvac 24 vacuum manifold and a vacuum was applied. (5) 700ul of Buffer RW1 was added to the column, and allowed to filter through the column. (6) 500ul of Buffer RPE was added onto the column and allowed to filter through. (7) The column was transferred to a 2ml collection tube, and centrifuged for 30 sec at 10,000 rpm. (8) 25µL of RNAase-free water was added to the membrane and the column was centrifuged for 30sec at 10,000rpm to elute the RNA. All centrifugation steps in the rapid protocol were performed on a model 10MVSS VWR centrifuge.

Extracted RNA was reverse transcribed, and RNA and cDNA were quantified as described previously.

For Taqman assays, Taqman Core Reagent Kit, Gene Amp 10X PCR Buffer 1, Amp Erase Uracil N-glycosidase, and AmpliTaq Gold DNA Polymerase were purchased from Applied Biosystems, Foster City, CA. All other reagents were from commercial sources described in example 1. Glycerol was purchased from Sigma Chemical Co (St. Louis, Mo.) and Tween™ 20 from Eastman Organic Chemicals (Rochester, NY)

Mammaglobin primers (SEQ ID NO. 3 and SEQ ID NO. 4) were synthesized by Invitrogen Corp.(Carlsbad, CA) and the mammaglobin Taqman probe (SEQ ID NO. 7) from Epoch Biosciences (San Diego, CA). PBGD primers (SEQ ID NO 8 and SEQ ID NO. 9) were synthesized by QIAGEN Operon (Alameda, CA), and the probe (SEQ ID NO. 23) by SYNTHEGEN, LLC (Houston, TX). B305D primers (SEQ ID NO. 11 and SEQ ID NO. 12) were synthesized at Invitrogen Corp and the probe (SEQ ID NO. 13) by Applied Biosystems, Inc. For all Taqman probes, carboxyfluorescein (FAM) and carboxytetramethylrhodamine TAMRA) were used as the dye and quencher pair.

SEQ ID NO. 3 CAAACGGATG AA ACTCTGAG CAATGTTGA

Docket No: CDS-5009

SEQ ID NO. 4 TCTGTGAGCC AAAGG TCTTG CAGA

SEQ ID NO. 7 6-FAM - TGTTTATGCA ATTAATATAT GACAGCAGTC
TTTGTG- TAMRA

5 SEQ ID NO. 8 CTGAGGCACC TGGAAGGAGG

SEQ ID NO. 9 CATCTTCATG CTGGGCAGGG

SEQ ID NO. 23 6- FAM-CCTGAGGCAC CTGGAAGGAG GCTGCAG TGT-
TAMRA

10 SEQ ID NO. 11 TCTGATAAAG GCCGTACAAT G

SEQ ID NO. 12 TCACGACTTG CTGTTTTTGC TC

SEQ ID NO. 13 6-FAM-ATCAAAAACA AGCATGGCCTA CACC- TAMRA

For the Taqman assays, a master mix was prepared by adding 10µL of 10X PCR
15 Buffer#1, 14 µL of 25mM MgCl₂, 8µL of 10mM dNTP's, 1µL of AmpErase UNG
(1U/µL), 16µL of glycerol, 1µL of 1%w/v Tween 20, 0.75µL of AmpliTaq Gold
(5U/µL), 6µL of each primer from a 5µM stock, 0.4µL of a 10µM stock of probe,
and water to a final volume of 96µL. Master mix (48 µL) and 2µL of cDNA from
each node sample were added to each optical reaction microtiter plate well. After
20 capping with optical caps, the plates were processed in an ABI Prism 7900 HT
Sequence Detection System (Applied Biosystems, Inc., Foster City, CA).
Commercial reagents including Taqman PCR Core Reagent Kit, Taqman DNA
Template Reagents, and Taqman B-actin Detection Reagents were purchased from
Applied Biosystems (Foster City, CA) and the assay was performed according to the
25 protocol recommended by the manufacturer. A threshold value of 0.02 was used for
analysis with the mammaglobin assay, 0.03 for the B305D assay, 0.08 for the B-
actin assay, and 0.1 for the PBDG assay. The average of triplicate determination by
Taqman assays are shown in Tables 3 and 4.

Docket No: CDS-5009

A comparison of the real-time measurement of the housekeeping genes B-actin and PBDG are shown in Table 1 in 15 H&E negative node samples as well as 16 H&E positive nodes and two control cDNA samples are shown in Table 1. Results of these experiments indicate good agreement between Ct values obtained with the rapid method of this invention and standard commercial methods (QIAGEN) confirming that RNA extracted by the rapid RNA extraction protocol is capable of being reverse transcribed and PCR amplified to a similar extent as RNA extracted based on a standard commercial protocol.

Table 2 compares gene expression results based on Ct value for the breast cancer markers mammaglobin and B305D with the same breast node samples evaluated in Table 1. To determine whether a correlation can be made between H&E positive and H&E negative breast node samples based on real-time PCR results, the lowest Ct value for each marker in the H&E negative node samples was identified for both mammaglobin and B305D. An arbitrary, but conservative, cut-off of 2.5Ct values less than this lowest value among H&E negative samples was applied to the H&E positive samples. Samples that are shaded in the Table 2 have higher expression for these two breast markers, in that Ct value is at least 2.5Ct values lower than a Ct value observed among H&E negative samples. There was a good correlation in expression of these markers using both the rapid method of this invention and the commercial manufacturer recommended protocol (QIAGEN).

Using mammaglobin as a marker, Ct values in two samples among the 15 H&E nodes samples did not strictly correlate between the rapid method of this invention and QIAGEN RNA extraction protocols. In one of these samples, GCLNC-24, the difference in Ct value was 0.4 which is within experimental error. The second sample showed a larger difference in Ct. larger differences in Ct value with both mammaglobin and B305D may be due to difficulties in spectrally quantifying RNA in these samples as well as to possible sampling problems due to the well established

Docket No: CDS-5009

non-uniform distribution of metastases in nodes. (Cserni. 1999. Metastases in axillary sentinel lymph nodes in breast cancer as detected by intensive histopathological work up. J. Clin Pathol, 52:0922-924.)

5

Docket No: CDS-5009

Table 1: Comparison of RNA extracted by the rapid method and Prior Art Method as measured by Taqman Assays for the Housekeeping Genes B-actin and PBGD

	Sample	Template	B-Actin Ct Values		PBGD Ct Values	
			Qiagen	Rapid	Qiagen	Rapid
H&E Negative	1	GCLNN-10	20.7	21.9	32.6	34.1
	2	GCLNN-11	20.5	21.0	33.7	32.7
	3	GCLNN-12	21.7	22.5	32.8	32.8
	4	GCLNN-14	21.2	19.1	32.1	31.0
	5	GCLNN-15	19.8	21.3	32.5	39.9
	6	GCLNN-19	22.4	23.0	29.9	39.1
	7	GCLNN-20	17.8	18.6	31.7	30.9
	8	GCLNN-23	21.4	19.2	33.1	31.4
	9	GCLNN-24	21.1	26.1	36.9	37.0
	10	GCLNN-25	19.1	21.4	32.7	33.1
	11	GCLNN-26	23.8	26.7	40.0	39.7
	12	CBLNN-274	20.1	19.8	31.5	30.7
	13	CBLNN-257	20.4	22.6	32.8	34.0
	14	CBLNN-258	23.4	27.6	33.1	35.4
	15	CBLNN-262	23.0	23.6	33.0	32.1
H&E Positive	16	GCLNC-1	18.7	18.5	28.8	29.1
	17	GCLNC-5	20.4	21.6	31.3	31.0
	18	GCLNC-6	20.0	19.9	31.4	30.2
	19	GCLNC-7	20.6	28.2	29.3	32.1
	20	GCLNC-11	18.9	19.5	31.2	30.9
	21	GCLNC-12	26.1	23.8	32.6	31.2
	22	GCLNC-13	20.5	18.9	30.8	27.7
	23	GCLNC-16	20.5	22.7	31.1	30.5
	24	GCLNC-17	22.1	19.4	40.0	29.1
	25	GCLNC-20	20.0	19.5	31.9	29.6
	26	GCLNC-21	21.1	21.3	33.5	32.2
	27	GCLNC-22	24.2	22.9	29.3	32.1
	28	GCLNC-23	18.4	20.7	29.6	35.5
	29	GCLNC-24	19.1	19.5	29.4	30.1
	31	ALNC-1	21.1	19.4	30.2	30.8
	32	ALNN-2	18.9	19.5	31.4	30.0
	33	CBLNN-247	18.2	20.7	31.5	31.1
Cut Off Values			33.7	37.5	36.1	36.1

Table 2: Comparison of RNA extracted by the rapid method and by the Prior Art Method as measured by by Taqman Assays for the Cancer Genes Mammaglobin and B305D

	Sample	Template	MG Ct Values		B305D-A/C Ct Values	
			Qiagen	Rapid	Qiagen	Rapid
H&E Negative	1	GCLNN-10	40.0	40.0	40.0	40.0
	2	GCLNN-11	40.0	40.0	40.0	40.0
	3	GCLNN-12	40.0	40.0	39.7	39.9
	4	GCLNN-14	40.0	40.0	39.6	39.9
	5	GCLNN-15	40.0	40.0	39.6	40.0
	6	GCLNN-19	40.0	40.0	40.0	40.0
	7	GCLNN-20	40.0	40.0	38.6	38.6
	8	GCLNN-23	40.0	40.0	39.0	40.0
	9	GCLNN-24	40.0	40.0	40.0	40.0
	10	GCLNN-25	40.0	40.0	40.0	40.0
	11	GCLNN-26	36.2	40.0	40.0	40.0
	12	CBLNN-274	40.0	40.0	39.3	39.5
	13	CBLNN-257	40.0	40.0	40.0	39.5
	14	CBLNN-258	40.0	40.0	40.0	40.0
	15	CBLNN-262	40.0	40.0	39.6	39.5
H&E Positive	16	GCLNC-1	20.8	22.8	27.1	27.9
	17	GCLNC-5	39.7	40.0	29.2	36.6
	18	GCLNC-6	23.7	22.7	26.8	26.3
	19	GCLNC-7	27.9	28.3	30.0	31.0
	20	GCLNC-11	32.5	32.8	28.2	28.4
	21	GCLNC-12	32.4	30.7	38.6	36.1
	22	GCLNC-13	40.0	39.7	38.9	38.3
	23	GCLNC-16	20.8	21.7	36.5	35.9
	24	GCLNC-17	30.7	28.7	32.3	24.9
	25	GCLNC-20	30.6	30.1	26.8	27.8
	26	GCLNC-21	40.0	38.9	34.0	35.2
	27	GCLNC-22	39.4	36.3	36.4	31.4
	28	GCLNC-23	33.0	33.8	30.1	30.8
	29	GCLNC-24	37.6	37.2	28.0	29.7
	31	ALNC-1	23.8	24.8	25.6	28.6
	32	ALNN-2	27.7	27.6	23.4	23.9
	33	CBLNN-247	35.9	38.2	29.8	30.1

Cut Off Values = 33.7 37.5 36.1 36.1

Docket No: CDS-5009

Example 2. Evaluation of a QIAshredder Column

It is the purpose of this example to demonstrate that the QIAshredder column can be used after homogenization of lymph node tissue, and that centrifugation of the homogenate for only 30 seconds yields acceptable RNA extraction.

5

A breast axillary node that was H&E -negative was obtained from Genomics Collaborative (Cambridge, MA). A 490mg slice of tissue was mixed with 5.5 ml of Buffer RLT, and the node was homogenized. Six hundred microliters of homogenate was removed and RNA was extracted using the rapid extraction described in

10 Example 1. As a control, another 600 μ L of homogenate was treated using the same rapid RNA extraction method except that a 3min centrifugation was used instead of the QIAshredder column. The RNA from each reaction was reverse transcribed and quantified as described in Example 1. Real-time PCR was performed using Taqman probes.

15

Results of these studies indicated similar RNA yields of 0.31 μ g/ μ L in the protocol involving the QIAshredder protocol as compared to 0.34 μ g/ μ L in the centrifugation protocol. Identical 260/280nm ratio's of 2.1 were obtained in both samples indicating high quality RNA. Real-time PCR assays also indicated similar Ct

20 values for Mammaglobin with RNA extracted by the QIAshredder column as compared to the protocol involving (22.5 versus 22.0), as well as with B305D (26.7 versus 26.6) as well as with the housekeeping genes PBDG (29.1 versus 29.0), and B-actin (SEQ ID NO 19, 18.5 versus 18.2).

25 In summary, this experiment indicates similar RNA yield, quality and real-time PCR amplification results with the protocol involving a QIAshredder instead of a 3min centrifugation step after homogenization. This reduces the time required to perform RNA extraction by 2.5 minutes, which is important in developing protocols suitable for intraoperative and other applications in which fast results are required to impact

Docket No: CDS-5009

patient care.

Example 3. Effect of Diluting of Lysate on RNA yield and precision

The following example illustrates the effect of dilution of the lysate on RNA yield and precision of RNA recovery. In this example, 20mg, 10mg or 5 mg of frozen pig node tissue was cut. Tissue was ground using a 50ml Disposable Tissue Grinder for 30 sec; and the samples were vortexed. Three replicates were performed for each node weight. One ml of each lysate was transferred to a 1.7 ml microcentrifuge tube which was then centrifuged in an Eppendorf microfuge at maximum speed (14,000RPM) for 30 sec. Seven hundred microliters was centrifuged through a QIAshredder column and the samples were centrifuged at maximum speed for 30 sec. The column was then washed, the column dried, and RNA eluted as described previously in Example 2, Part I steps 4-8. RNA was quantified by means of a Gene Spec II. Results of these studies are shown in Tables 3 and 4, and illustrates excellent precision at an initial node tissue weight of 2.5mg homogenized in 600ml of homogenization buffer, as compared to samples with higher weights including 10mg and 20mg node tissue. These results indicate that improved precision of RNA recovery can be obtained at lower node weight. Also, it would be expected that more dilute homogenates would be filtered faster through the column, and that there would be even fewer potential problems with filter clogging.

Docket No: CDS-5009

Table 3. Effect of Dilution of Lysate on RNA recovery





Table 4. Effect of Dilution.



Example 4: TaqMan testing of RNA from lymph node samples

Reverse transcription

cDNA was synthesized from 20 ug total RNA at 42 C for 60 min in a 250ul reaction
 5 containing 50 mM Tris-Cl pH 8.3, 75 mM KCl, 3 mM MgCl₂, 10 mM DTT, 2000U
 SuperscriptTM, 400U Rnasin, 2ug oligo dT primer, 0.08 mg/ml BSA and 0.2 mM each
 dACT, dCTP, dGTP and TTP. The resulting cDNA was diluted 1:8 in water.

PCR amplification

10 2ul of each diluted cDNA was amplified on the Applied Biosystems 7900 with a
 total volume per well of 50 ml. The amplification was carried out with the
 following protocol: 50°C for 2 min and 95°C for 10 min, then 40 cycles of 95°C for
 15 seconds (sec), 60°C for 1 min (58°C for B726) and 68°C for 1 min. The reaction
 mix components and final concentrations are as follows: 1X TaqMan buffer A
 15 (ABI), 3.5 mM MgCl₂, 0.2 mM of dGTP, dCTP and cATP, 0.4 mM dUTP (ABI),
 0.01 u/ul AmpErase UNG (ABI), 0.0375U/ul AmpliTaqTM Gold (ABI), 8% Glycerol
 (Sigma), 0.01% Tween 20 (KodakTM), 0.3 uM each primer (Invitrogen), RNase/DNase
 free water (Sigma) and 40 nM 5'-FAM, 3' TAMRA oligonucleotide probe
 (Synthegen). A threshold of signal relative to an internal reference dye (ROX) was
 20 assigned to each primer/probe set in the logarithmic linear range of the increase in
 fluorescent intensity caused by the increase in PCR product. The resulting Ct values
 were utilized to determine relative expression of the genes across all samples tested.

Primers utilized:

25

Mammaglobin

Primer 1, SEQ ID NO 3 - CAAACGGATGAAACTCTGAGCAATGTTGA

Primer 2, SEQ ID NO 4 - TCTGTGAGCCAAAGGTCTTGCAGA

Docket No: CDS-5009

Probe, SEQ ID NO 7 - FAM-

TGTTTATGCAATTAATATATGACAGCAGTCTTTGTG-TAMRA

B305D

5 Primer 1, SEQ ID NO 11 - TCTGATAAAGGCCGTACAATG

Primer 2, SEQ ID NO 12 - TCACGACTTGCTGTTTTTGCTC

Probe, SEQ ID NO 24 - FAM ATCAAAAAACAAGCATGGCCTCACACC

PBGD

10 Primer 1, SEQ ID NO 25 - CTGCTTCGCTGCATCGCTGAAA

Primer 2, SEQ ID NO 26 - CAGACTCCTCCAGTCAGGTACA

Probe, SEQ ID NO 23 - FAM-CCTGAGGCACCTGGAAGGAGGCTGCAGTGT-TAMRA

15 PIP

Primer 1, SEQ ID NO 27 - GCTTGGTGGTTAAACTTACC

Primer 2 SEQ ID NO 28, - TGAACAGTTCTGTTGGTGTA

Probe, SEQ ID NO 29 - FAM-CTGCCTGCCTATGTGACGACAATCCGG-TAMRA

20

GABA

Primer 1, SEQ ID NO 30 - CAATTTTGGTGGAGAACCCG

Primer 2, SEQ ID NO 31 - GCTGTCGGAGGTATATGGTG

Probe, SEQ ID NO 32 - FAM CATTTCAGAGAGTAACATGGACTACACA

25 TAMRA

B726

Primer 1, SEQ ID NO 33 - GCAAGTGCCAATGATCAGAGG

Primer 2, SEQ ID NO 34- ATATAGACTCAGGTATACACACT

Docket No: CDS-5009

Probe, SEQ ID NO 35 - FAM TCCCATCAGAATCCAAACAAGAGGAAG

Results:

- 5 The attached table contains the Ct values obtained for 11 histology-negative and 24 histology-positive lymph nodes. A cut-off for positivity for each marker was determined by subtracting 2.5 cycles from the lowest Ct observed with the histology-negative samples. Positive samples are denoted shading. Based on these criteria, detection rates for the individual markers were determined. Optimal
- 10 complementation is seen with mammaglobin combined with B305D A/C, a combination that detects 23/24 histology-positive samples. A combination of mammaglobin, B305D and GABA detected all 24 samples. See Table 5.

Docket No: CDS-5009

Sample	MG	PIP	B305D	GABA	B726	PBGD	Actin
N1	50.0	36.2	35.8	36.5	50.0	30.6	21.6
N2	50.0	38.7	39.0	37.0	50.0	30.4	21.1
N3	50.0	37.8	39.6	39.1	50.0	28.6	23.1
N4	36.3	36.4	39.8	40.0	50.0	32.4	19.8
N5	48.4	40.0	40.0	40.0	50.0	35.8	21.1
N6	46.9	34.8	38.9	36.1	50.0	28.5	20.3
N7	50.0	35.2	40.0	38.4	50.0	27.8	20.3
N8	38.8	38.5	36.9	37.5	50.0	29.4	21.7
N9	50.0	37.0	37.9	36.4	50.0	30.6	20.1
N10	50.0	40.0	38.1	40.0	50.0	31.0	22.3
N11	50.0	39.4	35.8	39.9	50.0	29.4	19.7
C1	25.7	25.1	25.9	36.7	25.8	28.8	21.7
C2	36.0	32.5	35.6	27.2	41.5	29.0	23.2
C3	26.2	30.7	26.8	37.2	35.2	29.9	24.1
C4	25.7	22.9	29.7	30.4	37.9	27.3	20.5
C5	25.9	23.3	34.8	38.0	39.1	27.8	21.1
C6	29.7	36.6	28.8	29.9	40.4	26.2	20.4
C7	35.0	33.9	25.4	23.4	42.7	27.1	22.5
C8	22.7	25.3	29.3	38.0	37.8	27.8	18.9
C9	19.2	26.6	27.7	36.3	36.0	26.2	21.1
C10	20.2	30.0	27.4	33.8	39.7	29.4	23.2
C11	30.4	33.2	29.0	26.2	40.7	27.9	21.4
C12	18.6	27.7	24.3	36.0	33.6	28.8	19.4
C13	23.4	17.0	23.1	38.2	24.3	27.4	21.8
C14	35.9	30.7	28.0	39.8	42.3	31.3	21.7
C15	27.0	28.1	25.6	34.6	25.5	29.7	21.1
C16	32.5	34.5	36.9	23.5	46.1	27.7	27.6
C17	38.8	35.7	27.5	26.3	37.4	28.8	29.7
C18	36.9	34.8	29.2	24.7	50.0	26.4	22.0
C19	31.6	22.9	24.5	38.0	39.1	28.9	24.7
C20	22.3	27.6	21.4	32.6	32.9	27.0	21.6
C21	19.3	24.8	25.1	36.2	28.9	28.4	21.3
C22	32.2	28.4	25.8	37.9	37.3	30.6	23.5
C23	20.7	24.2	27.2	37.4	27.4	28.9	22.6
C24	33.0	27.1	32.6	24.3	40.4	26.4	20.0
cutoff	33.8	32.3	33.3	33.6	37.5		

Positive samples detected

19/24

17/24

21/24

10/24

11/24

Table 5

5

10

Docket No: CDS-5009

Example 5 Multiple Markers

5 Reverse transcription

cDNA was synthesized from 20 ug total RNA at 42 C for 60 min in a 100ul reaction utilizing anchored oligo dT primers. 20 ng of each diluted cDNA was amplified on the Applied Biosystems 7900 with a total volume per well of 50 ml. The amplification was carried out with the following protocol: 94 C for 2 min then 50 cycles of: 94 C for 15 sec, 62 C for 15 sec and 72 C for 45 sec. The reaction mix components and final concentrations are as follows: 1X PCR Buffer (ABI), 2.5 mM final MgCl₂, 0.05 mM of dGTP, dCTP and cATP and TTP, 0.05U/ul Taq polymerase containing a 5-20X excess of 2 anti-Taq antibodies (TP4-9 and TP1-12), 0.2 uM each primer (Invitrogen), RNase/DNase free water (Sigma) and a 1:20000 dilution of a SYBR Green stock (Sigma). Amplification primer sequences were identical to those employed in TaqMan testing of lymph nodes. A common threshold of 275 fluorescence units was utilized to determine Ct values. Ct values were then converted into quantitative gene expression and plotted relative to a cut-off that was assumed to approximate the minimum level of gene expression required to differentiate metastatic lymph nodes from background expression .

]Primers utilized:

Mammaglobin

25 Primer 1, SEQ ID NO 3 - CAAACGGATGAACTCTGAGCAATGTTGA
Primer 2, SEQ ID NO 4 - TCTGTGAGCCAAAGGTCTTGCAGA

B305D

Primer 1, SEQ ID NO 11 - TCTGATAAAGGCCGTACAATG
30 Primer 2, SEQ ID NO 12 - TCACGACTTGCTGTTTTTGCTC

Docket No: CDS-5009

This data displayed in Table 6-8 shows the power of combining genes
mammaglobin and B305D signatures for obtaining coverage of the vast majority of
breast cancer samples.


5

Docket No: CDS-5009

Table 6.

5

Expression of Mammaglobin in Different Breast Cell Types



10


15

Docket No: CDS-5009

Table 7

5

Expression of B305D in Different Breast Cell Types




10

Docket No: CDS-5009

Table 8

Combined Expression of Mammaglobin & B305D in Different Breast Cell Types



Docket No: CDS-5009

Example 5: Rapid RT-PCR from purified RNA

- 5 Mammaglobin in vitro transcript RNA (100,000 copies) was diluted into 20 ng white blood cell RNA and subjected to one-step RT-PCR amplification in a Cepheid Smart Cycler II in a reaction containing the following: 50 mM Bicine / KOH, pH 8.2, 115 mM Potassium Acetate, 8% v/v Glycerol, 0.2 mg/ml BSA, 150 mM Trehalose, 0.2% v/v Tween 20, 0.2 mM Tris-Cl pH 8, 3.5 mM MnSO₄, 0.3 mM
10 dATP, 0.3 mM dCTP, 0.3 mM dGTP, 0.3 mM TTP, 100 U/ml Tth, 20 ug/ml Antibody TP6-25.3, 0.45 uM Mammaglobin Scorpion and 0.5 uM Mammaglobin Primer.

- The PCR conditions were as follows (the asterisk denoted the portion of cycling in
15 which the fluorescence was measured:

Condition A - 95 C for 3 sec, 60 C for 7 min, 70 C for 2 min, then 40 cycles of: 95 C for 3 sec, 54 C for 6 sec* and 68 C for 6 sec – total time = 34 min

- 20 Condition B - 95 C for 3 sec, 60 C for 3 min, then 40 cycles of: 95 C for 3 sec, 54 C for 6 sec* and 68 C for 6 sec – total time = 29 min

Condition C - 95 C for 3 sec, 60 C for 3 min, then 40 cycles of: 95 C for 3 sec, 58 C for 6 sec* and 68 C for 6 sec – total time = 26 min

25

The results of this assay are shown in the following table:

Docket No: CDS-5009

Cond	template	mean Ct	Total Assay Time (min)
A	IVT-RNA	30.47	34.00
A	NT	50.00	"
B	IVT-RNA	30.26	29.00
B	NT	50.00	"
C	IVT-RNA	31.15	26.00
C	NT	50.00	"

As can be seen from the data, reducing the length of the reverse transcription step from 9 min to 3 min has no discernable impact on assay performance. Increasing the annealing temperature to 58 C has only a minimal impact of amplification and reduces the total RT-PCR time to 26 minutes. Further reductions in cycling time can be accomplished by primer and probe construction with higher optimal assay temperatures.

10

Example 6: Rapid multiplex RT-PCR assay for detection of lymph node metastasis (Prophetic)

In this example, purified lymph node RNA (≥ 100 ng total RNA per 100ul reaction utilized to reduce the number of cycles required to generate a positive result) is subjected to one-step RT-PCR amplification in a Cepheid Smart Cycler II. The reaction contains the following components: 50 mM Bicine / KOH, pH 8.2, 115 mM Potassium Acetate, 8% v/v Glycerol, 0.2 mg/ml BSA, 150 mM Trehalose, 0.2% v/v Tween 20, 0.2 mM Tris-Cl pH 8, 3.5 mM MnSO₄ or Mn(acetate), 0.3 mM dATP, 0.3 mM dCTP, 0.3 mM dGTP, 0.3 mM TTP, 100 U/ml Tth, 20 ug/ml Antibody TP6-25.3, 0.4-0.8 uM each Scorpion and primer. Scorpion and primer sequences are designed such that non-specific priming events are minimized and such that they are

Docket No: CDS-5009

compatible with multiplexed amplification / detection with the following rapid RT-PCR profile:

95 C for 3 sec, 60 C for 3 min, then 32 cycles of: 95 C for 1 sec, 65 C for 6 sec* and 72 C for 3 sec – total amplification time = 16 min.

5

The rapid sample preparation method described is utilized that takes ≤ 4 min to complete (refer to sample prep example in application), leading to a total assay time of ≤ 20 min.

- 10 For detection of metastasis of breast cancer to lymph nodes, the following genes are amplified: mammaglobin, B305D and PBGD. Mammaglobin and/or B305D signal is indicative of metastasis, while PBGD is utilized as a housekeeping gene. In this example, Ct values for cancer-positive lymph nodes range from 12-30 for mammaglobin and 12-28 for B305D – detection of either gene in this Ct range is
- 15 indicative of a cancer-positive lymph node. In this example, Ct values for the control gene PBGD must fall in the range of 24-30 cycles in order for a cancer-negative result to be considered valid.

- For detection of metastasis of melanoma cancer to lymph nodes, the same conditions
- 20 are utilized, with the following exceptions: (1) lymph nodes suspected of containing melanoma are processed and (2) the following genes are amplified: MART1, tyrosinase and PBGD. MART1 and/or tyrosinase signal is indicative of metastasis of melanoma to the lymph nodes, while PBGD is utilized as a housekeeping gene. In this example, Ct values for cancer-positive lymph nodes range from 18-30 for
- 25 tyrosinase and 12-27 for MART1 – detection of either gene in this range is indicative of a cancer-positive lymph node. In this example, Ct values for the control gene PBGD must fall in the range of 24-30 cycles in order for a cancer-negative result to be considered valid.

SEQUENCE LISTING

<110> Veridex, LLC

<120> INTRAOPERATIVE LYMPH NODE ASSAY

<130> 1011-4416CA

<160> 35

<170> PatentIn version 3.1

<210> 1

<211> 23

<212> DNA

<213> Artificial

<220>

<223> Artificial sequence is a Primer

<400> 1

atgactgcct tgcctcctca gta

23.

<210> 2

<211> 26

<212> DNA

<213> Artificial

<220>

<223> Artificial sequence is a Primer

<400> 2

ggctgttgct tggacttctc taaaga

26

<210> 3

<211> 29

<212> DNA

<213> Artificial

<220>

<223> Artificial sequence is a Primer

<400> 3

caaacggatg aaactctgag caatgttga

29

<210> 4

<211> 24

<212> DNA

<213> Artificial

<220>

<223> Artificial sequence is a Primer

<400> 4

tctgtgagcc aaaggtcttg caga

24

<210> 5
 <211> 23
 <212> DNA
 <213> Artificial

<220>
 <223> Artificial sequence is a Primer

<400> 5
 ggccaacaaa gctcaggaca aca

23

<210> 6
 <211> 25
 <212> DNA
 <213> Artificial

<220>
 <223> Artificial sequence is a Primer

<400> 6
 gcagtgactt cgtcatttgg acgta

25

<210> 7
 <211> 36
 <212> DNA
 <213> Artificial

<220>
 <223> Artificial sequence is a Probe

<400> 7
 tgtttatgca attaatatat gacagcagtc tttgtg

36

<210> 8
 <211> 20
 <212> DNA
 <213> Artificial

<220>
 <223> Artificial sequence is a Primer

<400> 8
 ctgaggcacc tggaaggagg

20

<210> 9
 <211> 20
 <212> DNA
 <213> Artificial

<220>

<223> Artificial sequence is a Primer

<400> 9
catcttcatg ctgggcaggg

20

<210> 10
<211> 29
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 10
ctgaggcacc tggaaggagg ctgcagtgt

29

<210> 11
<211> 21
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 11
tctgataaag gccgtacaat g

21

<210> 12
<211> 22
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 12
tcacgacttg ctgtttttgc tc

22

<210> 13
<211> 26
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Probe

<400> 13
atcaaaaaac aagcatggcc tacacc

26

<210> 14
<211> 1155
<212> DNA

<213> Human

<400> 14

atggtggttg aggttgattc catgccggct gcctcttctg tgaagaagcc atttgggtctc 60
 aggagcaaga tgggcaagtg gtgctgccgt tgcttcccct gctgcaggga gagcggcaag 120
 agcaacgtgg gcacttctgg agaccaagac gactctgcta tgaagacact caggagcaag 180
 atgggcaagt ggtgccgcca ctgcttcccc tgctgcaggg ggagtggcaa gagcaacgtg 240
 ggcgcttctg gagaccacga cgactctgct atgaagacac tcaggaacaa gatgggcaag 300
 tgggtgctgcc actgcttccc ctgctgcagg gggagcagca agagcaaggt gggcgcttgg 360
 ggagactacg atgacagtgc cttcatggag cccagggtacc acgtccgtgg agaagatctg 420
 gacaagctcc acagagctgc ctggtggggg aaagtcccca gaaaggatct catcgtcatg 480
 ctcagggaca ctgacgtgaa caagcaggac aagcaaaaga ggactgctct acatctggcc 540
 tctgccaatg ggaattcaga agtagtaaaa ctctgctgg acagacgatg tcaacttaat 600
 gtccttgaca acaaaaagag gacagctctg ataaaggccg tacaatgcca ggaagatgaa 660
 tgtgcgttaa tggtgctgga acatggcact gatccaaata ttccagatga gtatggaaat 720
 accactctgc actacgtat ctataatgaa gataaattaa tggccaaagc actgctctta 780
 tatggtgctg atatcgaatc aaaaaacaag catggcctca caccactgtt acttgggtga 840
 catgagcaaa aacagcaagt cgtgaaatth ttaattaaga aaaaagcgaa tttaaagtga 900
 ctggatagat atggaaggac tgctctcata cttgctgtat gttgtggatc agcaagtata 960
 gtcagccttc tacttgagca aaatattgat gtatcttctc aagatctatc tggacagacg 1020
 gccagagagt atgctgtttc tagtcatcat catgtaatth gccagttact ttctgactac 1080
 aaagaaaaac agatgctaaa aatctcttct gaaaacagca atccagaaaa tgtctcaaga 1140
 accagaaata aataa 1155

<210> 15

<211> 3865

<212> DNA

<213> Human

<400> 15

tccgagctga ttacagacac caaggaagat gctgtaaaga gtcagcagcc acagccctgg 60
 ctagctggcc ctgtgggcat ttattagtaa agttttaatg acaaaagctt tgagtcaaca 120
 caccgtggg taattaacct ggtcatcccc accctggaga gccatcctgc ccatgggtga 180
 tcaaagaagg aacatctgca ggaacacctg atgaggctgc acccttggcg gaaagaacac 240

ctgacacagc tgaaagcttg gtggaaaaaa cacctgatga ggctgcaccc ttggtggaaa	300
gaacacctga cacggctgaa agcttggtgg aaaaaacacc tgatgaggct gcaccccttg	360
tggaggggaa atctgacaaa attcaatgtt tggagaaagc gacatctgga aagttcgaac	420
agtcagcaga agaaacacct agggaaatta cgagtcctgc aaaagaaaca tctgagaaat	480
ttacgtggcc agcaaaagga agacctagga agatcgcatg ggagaaaaaa gaagacacac	540
ctaggggaaat tatgagtcct gcaaaagaaa catctgagaa atttacgtgg gcagcaaaag	600
gaagacctag gaagatcgca tgggagaaaa aagaaacacc tgtaaagact ggatgcgtgg	660
caagagtaac atctaataaa actaaagttt tggaaaaagg aagatctaag atgattgcat	720
gtcctacaaa agaatcatct acaaaagcaa gtgccaatga tcagagggtc ccatcagaat	780
ccaaacaaga ggaagatgaa gaatattctt gtgattctcg gagtctcttt gagagtcttg	840
caaagattca agtgtgtata cctgagtcta tatatcaaaa agtaatggag ataaatagag	900
aagtagaaga gcctcctaag aagccatctg ccttcaagcc tgccattgaa atgcaaaaact	960
ctgttccaaa taaagccttt gaattgaaga atgaacaaac attgagagca gatccgatgt	1020
tcccaccaga atccaaacaa aaggactatg aagaaaattc ttgggattct gagagtctct	1080
gtgagactgt ttcacagaag gatgtgtgtt taccgaaggc tacacatcaa aaagaaatag	1140
ataaaataaa tggaaaatta gaagagtctc ctaataaaga tgggtcttctg aaggctacct	1200
gcggaatgaa agtttctatt ccaactaaag ccttagaatt gaaggacatg caaactttca	1260
aagcagagcc tccggggaag ccatctgcct tcgagcctgc cactgaaatg caaaagtctg	1320
tcccaaataa agccttgga tggaaaaatg aacaaacatt gagagcagat gagatactcc	1380
catcagaatc caaacaaaag gactatgaag aaagttcttg ggattctgag agtctctgtg	1440
agactgtttc acagaaggat gtgtgtttac ccaaggctrc rcatcaaaaa gaaatagata	1500
aaataaatgg aaaattagaa gggctctcctg ttaaagatgg tcttctgaag gctaactgcg	1560
gaatgaaagt ttctattcca actaaagcct tagaattgat ggacatgcaa actttcaaag	1620
cagagcctcc cgagaagcca tctgccttcg agcctgccat tgaaatgcaa aagtctgttc	1680
caaataaagc cttggaattg aagaatgaac aaacattgag agcagatgag atactcccat	1740
cagaatccaa acaaaaggac tatgaagaaa gttcttggga ttctgagagt ctctgtgaga	1800
ctgtttcaca gaaggatgtg tggtttacca aggctrrcca tcaaaaagaa atagataaaa	1860
taaatggaaa attagaagag tctcctgata atgatggttt tctgaaggct cctgcagaa	1920

tgaaagtttc tattccaact aaagccttag aattgatgga catgcaaact ttcaaagcag	1980
agcctcccga gaagccatct gccttcgagc ctgcoattga aatgcaaaag tctgttccaa	2040
ataaagcctt ggaattgaag aatgaacaaa cattgagagc agatcagatg ttcccttcag	2100
aatcaaaaaca aaagaasgtt gaagaaaatt cttgggattc tgagagtctc cgtgagactg	2160
tttcacagaa ggatgtgtgt gtacccaagg ctacacatca aaaagaaatg gataaaataa	2220
gtggaaaatt agaagattca actagcctat caaaaatctt ggatacagtt cattcttgtg	2280
aaagagcaag ggaacttcaa aaagatcact gtgaacaacg tacaggaaaa atggaacaaa	2340
tgaaaaagaa gttttgtgta ctgaaaaaga aactgtcaga agcaaaagaa ataaaatcac	2400
agttagagaa caaaaaagtt aaatgggaac aagagctctg cagtgtgaga ttgactttaa	2460
accaagaaga agagaagaga agaaatgccg atatattaaa tgaaaaaatt agggaagaat	2520
taggaagaat cgaagagcag cataggaaag agttagaagt gaaacaacaa cttgaacagg	2580
ctctcagaat acaagatata gaattgaaga gtgtagaaag taatttgaat caggtttctc	2640
acactcatga aaatgaaaat tatctcttac atgaaaattg catgttgaaa aaggaaattg	2700
ccatgctaaa actggaaata gccacactga aacaccaata ccaggaaaag gaaaataaat	2760
actttgagga cattaagatt ttaaaagaaa agaatgctga acttcagatg accctaaaac	2820
tgaaagagga atcattaact aaaagggcat ctcaatatag tgggcagctt aaagttctga	2880
tagctgagaa cacaatgctc acttctaaat tgaaggaaaa acaagacaaa gaaatactag	2940
aggcagaaat tgaatcacac catcctagac tggcttctgc tgtacaagac catgatcaaa	3000
ttgtgacatc aagaaaaagt caagaacctg ctttccacat tgcaggagat gcttgtttgc	3060
aaagaaaaat gaatgttgat gtgagtagta cgatatataa caatgaggtg ctccatcaac	3120
cactttctga agctcaaagg aaatccaaaa gcctaaaaat taatctcaat tatgcmggag	3180
atgctctaag agaaaataca ttggtttccag aacatgcaca aagagaccaa cgtgaaacac	3240
agtgtcaaat gaaggaagct gaacacatgt atcaaaacga acaagataat gtgaacaaac	3300
acactgaaca gcaggagtct ctagatcaga aattatttca actacaaagc aaaaatatgt	3360
ggcttcaaca gcaattagtt catgcacata agaaagctga caacaaaagc aagataacaa	3420
ttgatattca ttttcttgag aggaaaatgc aacatcatct cctaaaagag aaaaatgagg	3480
agatatttaa ttacaataac catttaaaaa accgtatata tcaatatgaa aaagagaaag	3540
cagaaacaga aaactcatga gagacaagca gtaagaaact tcttttgag aaacaacaga	3600
ccagatcttt actcacaact catgctagga ggccagtcct agcatcacct tatgttgaaa	3660

atcttaccaa tagtctgtgt caacagaata cttatttttag aagaaaaatt catgatttct 3720
 tcctgaagcc tacagacata aaataacagt gtgaagaatt acttggttcac gaattgcata 3780
 aagctgcaca ggattcccat ctaccctgat gatgcagcag acatcattca atccaaccag 3840
 aatctcgctc tgtcactcag gctgg 3865

<210> 16
 <211> 3282
 <212> DNA
 <213> Human

<400> 16
 gggacagggc tgaggatgag gagaaccctg gggacccaga agaccgtgcc ttgcccggaa 60
 gtcttgctg taggcctgaa ggacttgccc taacagagcc tcaacaacta cctgggtgatt 120
 cctacttcag ccccttggtg tgagcagctt ctcaacatga actacagcct ccacttggcc 180
 ttcgtgtgtc tgagtctctt cactgagagg atgtgcatcc aggggagtca gttcaacgtc 240
 gaggtcggca gaagtgacaa gctttccctg cctggctttg agaacctcac agcaggatat 300
 aacaaatttc tcaggcccaa ttttggtgga gaaccgtac agatagcgtt gactctggac 360
 attgcaagta tctctagcat ttcagagagt aacatggact acacagccac catatacctc 420
 cgacagcgtt ggatggacca gcggtggtg tttgaaggca acaagagctt cactctggat 480
 gccgcctcg tggagttcct ctgggtgcca gatacttaca ttgtggagtc caagaagtcc 540
 ttctccatg aagtcactgt gggaaacagg ctcatccgcc tcttctccaa tggcacggtc 600
 ctgtatgcc tcagaatcac gacaactgtt gcatgtaaca tggatctgtc taaatacccc 660
 atggacacac agacatgcaa gttgcagctg gaaagctggg gctatgatgg aaatgatgtg 720
 gaggttcacct ggctgagagg gaacgactct gtgcgtggac tggaacacct gcggcttgct 780
 cagtacacca tagagcggta tttcacctta gtcaccagat cgcagcagga gacaggaaat 840
 tacactagat tggctttaca gtttgagctt cggaggaatg ttctgtatct cattttggaa 900
 acctacgttc ctteccacttt cctgggtggtg ttgtcctggg tttcattttg gatctctctc 960
 gattcagtc ctgcaagaac ctgcattgga gtgacgaccg tggttatcaat gaccacactg 1020
 atgatcgggt cccgcacttc ttttcccaac accaactgct tcatcaaggc catcgatgtg 1080
 tacctgggga tctgcttttag ctttgtgttt ggggccttgc tagaatatgc agttgctcac 1140
 tacagttcct tacagcagat ggcagccaaa gataggggga caacaaagga agtagaagaa 1200
 gtcagtatta ctaatatcat caacagctcc atctccagct ttaaacggaa gatcagcttt 1260

gccagcattg aaatttccag cgacaacggt gactacagtg acttgacaat gaaaaccagc 1320
 gacaagttca agtttgtctt ccgagaaaag atgggcagga ttgttgatta tttcacaatt 1380
 caaaacccca gtaatgttga tcactattcc aaactactgt ttcctttgat ttttatgcta 1440
 gccaatgtat tttactgggc atactacatg tttttttgag tcaatgttaa atttcttgca 1500
 tgccataggt cttcaacagg acaagataat gatgtaaag gtattttagg ccaagtgtgc 1560
 acccacatcc aatggtgcta caagtgactg aaataatatt tgagtctttc tgctcaaaga 1620
 atgaagctcc aaccattggt ctaagctgtg tagaagtcct agcattatag gatcttgtaa 1680
 tagaaacatc agtccattcc tttttcatct taatcaagga cattcccatg gagcccaaga 1740
 ttacaaatgt actcaggggt gtttattcgg tggctccctg gtttgcatth acctcatata 1800
 aagaatggga aggagaccat tgggtaacc tcaagtgtca gaagttgttt ctaaagtaac 1860
 tatacatggt ttttactaaa tctctgcagt gcttataaaa tacattgttg cctatttagg 1920
 gagtaacatt ttctagtttt tgtttctggt taaaatgaaa tatgggctta tgtcaattca 1980
 ttggaagtca atgcactaac tcaataccaa gatgagtttt taaataatga atattattta 2040
 ataccacaac agaattatcc ccaatttcca ataagtccta tcattgaaaa ttcaaata 2100
 agtgaagaaa aaattagtag atcaacaatc taaacaaatc cctcggttct aagatacaat 2160
 ggattcccca tactggaagg actctgaggc tttattcccc cactatgcat atcttatcat 2220
 tttattatta tacacacatc catcctaaac tatactaaag cccttttccc atgcatggat 2280
 ggaaatggaa gatttttttg taacttggtc tagaagtctt aatatgggct gttgccatga 2340
 aggcttgagc aattgagtc cttttctagc tgcctttatt cacatagtga tggggacta 2400
 aaagtactgg gttgactcag agagtcgctg tcattctgtc attgctgcta ctctaacact 2460
 gagcaacact ctcccagtg cagatccct gtatcattcc aagaggagca ttcattccct 2520
 tgctctaag atcaggaatg atgcttatta gaaaacaaac tgcttgacct aggaacaagt 2580
 ggcttagctt aagtaaactt ggctttgctc agatccctga tccttccagc tggctgctc 2640
 tgagtggctt atcccgcag agcaggagcg tgctggcctt gactactgaa ctttctgagt 2700
 aacaatgaga caggttacag aacctatgtt caggttgagg gtgagctgcc ctctccaaat 2760
 ccagccagag atgcacattc ctgggccagt ctgagccaac agtaccacaaa gtgatttttg 2820
 agtgtgccag ggtaaaggct tccagttcag cctcagttat tttagacaat ctgcctatct 2880
 ttaatttctt agcttcctgt tctaataaat gcacggcttt acctttcctg tcagaaataa 2940

accaaggctc taaaagatga tttcccttct gtaactccct agagccacag gttctcatte 3000
 cttttcccat tatacttctc acaattcagt ttctatgagt ttgatcacct gattttttta 3060
 acaaaatatt tctaacggga atgggtggga gtgctggtga aaagagatga aatgtggttg 3120
 tatgagccaa tcatatttgt gattttttta aaaaagttaa aaaggaaata tctgttctga 3180
 aacccactt aagcattggt tttatataaa aacaatgata aagatgtgaa ctgtgaaata 3240
 aatataccat attagctacc caccaaaaaa aaaaaaaaaa aa 3282

<210> 17
 <211> 503
 <212> DNA
 <213> Human

<400> 17
 gacagcggct tccttgatcc ttgccaccgc cgactgaaca ccgacagcag cagcctcacc 60
 atgaagttgc tgatggctct catgctggcg gccctctccc agcactgcta cgcaggctct 120
 ggctgcccct tattggagaa tgtgatttcc aagacaatca atccacaagt gtctaagact 180
 gaatacaaag aacttcttca agagttcata gacgacaatg ccactacaaa tgccatagat 240
 gaattgaagg aatgttttct taaccaaacg gatgaaactc tgagcaatgt tgaggtgttt 300
 atgcaattaa tatatgacag cagtctttgt gatttatttt aactttctgc aagacctttg 360
 gctcacagaa ctgcagggtg tgggtgagaaa ccaactacgg attgctgcaa accacacctt 420
 ctctttctta tgtcttttta ctacaaacta caagacaatt gttgaaacct gctatacatg 480
 tttattttta taaattgatg gca 503

<210> 18
 <211> 503
 <212> DNA
 <213> Human

<400> 18
 gacagcggct tccttgatcc ttgccaccgc cgactgaaca ccgacagcag cagcctcacc 60
 atgaagttgc tgatggctct catgctggcg gccctctccc agcactgcta cgcaggctct 120
 ggctgcccct tattggagaa tgtgatttcc aagacaatca atccacaagt gtctaagact 180
 gaatacaaag aacttcttca agagttcata gacgacaatg ccactacaaa tgccatagat 240
 gaattgaagg aatgttttct taaccaaacg gatgaaactc tgagcaatgt tgaggtgttt 300
 atgcaattaa tatatgacag cagtctttgt gatttatttt aactttctgc aagacctttg 360
 gctcacagaa ctgcagggtg tgggtgagaaa ccaactacgg attgctgcaa accacacctt 420

ctctttctta tgtcttttta ctacaaacta caagacaatt gttgaaacct gctatacatg 480

tttatttttaa taaattgatg gca 503

<210> 19
 <211> 1793
 <212> DNA
 <213> Human

<400> 19
 cgcgctccgcc ccgcgagcac agagcctcgc ctttgccgat ccgcgcgccg tccacacccg 60
 ccgcoagctc accatggatg atgatatcgc cgcgctcgtc gtcgacaacg gctccggcat 120
 gtgcaaggcc ggcttcgcgg gcgacgatgc cccccgggcc gtcttcccct ccatcgtggg 180
 gcgccccagg caccagggcg tgatgggtgg catgggtcag aaggattcct atgtgggcga 240
 cgaggccccag agcaagagag gcatactcac cctgaagtac cccatcgagc acggcatcgt 300
 caccaactgg gacgacatgg agaaaatctg gcaccacacc ttctacaatg agctgcgtgt 360
 ggctcccgag gagcaccgcc tgctgctgac cgaggcccc ctgaaccca aggccaaccg 420
 cgagaagatg acccagatca tgtttgagac cttcaacacc ccagccatgt acgttgctat 480
 ccaggctgtg ctatccctgt acgcctctgg ccgtaccact ggcatcgtga tggactccgg 540
 tgacgggggtc acccactg tgcccatcta cgaggggtat gccctcccc atgccatcct 600
 gcgtctggac ctggctggcc gggacctgac tgactacctc atgaagatcc tcaccgagcg 660
 cggctacagc ttcaccacca cggccgagcg ggaaatcgtg cgtgacatta aggagaagct 720
 gtgctacgtc gccctggact tcgagcaaga gatggccacg gctgcttcca gctcctccct 780
 ggagaagagc tacgagctgc ctgacggcca ggtcatcacc attggcaatg agcggttccg 840
 ctgccctgag gcactcttcc agccttcctt cctgggcatg gagtcctgtg gcatccacga 900
 aactaccttc aactccatca tgaagtgtga cgtggacatc cgcaaagacc tgtacgccaa 960
 cacagtgtg tctggcgcca ccaccatgta ccctggcatt gccgacagga tgcagaagga 1020
 gatcactgcc ctggcaccca gcacaatgaa gatcaagatc attgctcctc ctgagcgcaa 1080
 gtactcogtg tggatcggcg gctccatcct ggctcogctg tccacottcc agcagatgtg 1140
 gatcagcaag caggagtatg acgagtcggg cccctccatc gtccaccgca aatgcttcta 1200
 ggcggaactat gacttagttg cgttacacco tttcttgaca aaacctact tgcgcagaaa 1260
 acaagatgag attggcatgg ctttatttgt tttttttgtt ttgttttggt tttttttttt 1320
 tttttggctt gactcaggat ttaaaaactg gaacggtgaa ggtgacagca gtcggttgga 1380

gcgagcatcc cccaaagttc acaatgtggc cgaggacttt gattgcacat tggtgttttt 1440
 ttaatagtca ttccaaatat gagatgcatt gttacaggaa gtcccttgcc atcctaaaag 1500
 ccacccact tctctctaag gagaatggcc cagtcctctc ccaagtccac acaggggagg 1560
 tgatagcatt gctttcgtgt aaattatgta atgcaaaatt tttttaatct tcgccttaat 1620
 acttttttat tttgttttat tttgaatgat gagccttcgt gccccccctt cccctttttt 1680
 gtcccccaac ttgagatgta tgaaggcttt tgggtctcct gggagtgggt ggaggcagcc 1740
 agggcttacc tgtacactga cttgagacca gttgaataaa agtgcacacc tta 1793

<210> 20
 <211> 1377
 <212> DNA
 <213> Human

<400> 20
 cacacagcct actttccaag cggagccatg tctggtaacg gcaatgcggc tgcaacggcg 60
 gaagaaaaca gcccaaagat gagagtgatt cgcgtgggta cccgcaagag ccagcttgct 120
 cgcatacaga cggacagtgt ggtggcaaca ttgaaagcct cgtaccctgg cctgcagttt 180
 gaaatcattg ctatgtccac cacaggggac aagattcttg atactgcact ctctaagatt 240
 ggagagaaaa gcctgtttac caaggagctt gaacatgccc tggagaagaa tgaagtggac 300
 ctggttgttc actccttgaa ggacctgccc actgtgcttc ctctgggctt caccatcgga 360
 gccatctgca agcgggaaaa cctcatgat gctgttgtct ttcacccaaa atttggtggg 420
 aagaccctag aaaccctgcc agagaagagt gtggtgggaa ccagctccct gcgaagagca 480
 gccagctgc agagaaagtt cccgcatctg gagttcagga gtattcgggg aaacctcaac 540
 acccggttc ggaagctgga cgagcagcag gagttcagtg ccatcatcct agcaacagct 600
 ggctgcagc gcatgggctg gcacaaccgg gtggggcaga tcctgcaccc tgagaaatgc 660
 atgtatgctg tgggccaggg ggcttgggc gtggaagtgc gagccaagga ccaggacatc 720
 ttggatctgg tgggtgtgct gcacgatccc gagactctgc ttcgctgcat cgctgaaagg 780
 gccttcctga ggcacctgga aggaggctgc agtgtgccag tagcogtgca tacagctatg 840
 aaggatgggc aactgtacct gactggagga gtctggagtc tagacggctc agatagcata 900
 caagagacca tgcaggctac catccatgtc cctgcccagc atgaagatgg ccctgaggat 960
 gaccacagt tggtaggcct cactgctcgt aacattccac gagggcccca gttggctgcc 1020
 cagaacttgg gcatcagcct ggccaacttg ttgctgagca aaggagccaa aaacatcctg 1080

gatgttgac ggcagcttaa cgatgcccac taactgggtt gtggggcaca gatgcctggg 1140
 ttgctgctgt ccagtgccta catcccgggc ctcagtgcc cattctcact gctatctggg 1200
 gagtgattac ccogggagac tgaactgcag ggttcaagcc ttccagggat ttgcctcacc 1260
 ttggggcctt gatgactgcc ttgcctcttc agtatgtggg ggcttcatct ctttagagaa 1320
 gtccaagcaa cagcctttga atgtaaccaa tcctactaat aaaccagttc tgaaggt 1377

<210> 21
 <211> 1524
 <212> DNA
 <213> Human

<400> 21
 agcagacaga ggactctcat taaggaaggt gtctgtgcc ctgaccctac aagatgccaa 60
 gagaagatgc tcaattcatc tatggttacc ccaagaaggg gcacggccac tcttacacca 120
 cggctgaaga ggccgctggg atcgccatcc tgacagtgat cctgggagtc ttactgctca 180
 tcggctgttg gtattgtaga agacgaaatg gatacagagc cttgatggat aaaagtcttc 240
 atgttggcac tcaatgtgcc ttaacaagaa gatgcccaca agaagggtt gatcatcggg 300
 acagcaaagt gtctcttcaa gagaaaaact gtgaacctgt gggtcccaat gctccacctg 360
 cttatgagaa actctctgca gaacagtcac caccacctta ttcaccttaa gagccagcga 420
 gacacctgag acatgctgaa attatttctc tcacactttt gcttgaattt aatacagaca 480
 tctaattgtc tcctttggaa tgggtgtagga aaaatgcaag ccatctctaa taataagtca 540
 gtgttaaaat tttagtaggt ccgctagcag tactaatcat gtgaggaaat gatgagaaat 600
 attaaattgg gaaaactcca tcaataaatg ttgcaatgca tgatactatc tgtgccagag 660
 gtaatgttag taaatccatg gtgttatatt ctgagagaca gaattcaagt gggattcttg 720
 gggccatcca atttctcttt acttgaaatt tggctaataa caaactagtc aggttttcga 780
 accttgaccg acatgaactg tacacagaat tgttccagta ctatggagtg ctcacaaagg 840
 atacttttac aggttaagac aaagggttga ctggcctatt tatctgatca agaacatgtc 900
 agcaatgtct ctttgtgctc taaaattcta ttatactaca ataatatatt gtaaagatcc 960
 tatagctctt tttttttgag atggagtttc gcttttggtg ccagggctgg agtgcaatgg 1020
 cgcgatcttg gctcaccata acctccgcct ccaggttca agcaattctc ctgccttagc 1080
 ctctgagta gctgggatta caggcgtgcg ccactatgcc tgactaattt tgtagtttta 1140
 gtagagacgg ggtttctcca tgttggtcag gctggtctca aactcctgac ctcagggtgat 1200

ctgcccgcct cagcctccca aagtgcctgga attacaggcg tgagccacca cgccctggctg 1260
gatcctatat cttaggtaag acatataacg cagtctaatt acattttcact tcaaggctca 1320
atgctattct aactaatgac aagtattttc tactaaacca gaaattggta gaaggattta 1380
aataagtaaa agctactatg tactgcctta gtgctgatgc ctgtgtactg ccttaaattgt 1440
acctatggca atttagctct cttgggttcc caaatccctc tcacaagaat gtgcagaaga 1500
aatcataaag gatcagagat tctg 1524

<210> 22
<211> 2384
<212> DNA
<213> Human

<400> 22
tattgagttc ttcaaacatt gtagcctctt tatggtctct gagaaataac taccttaaac 60
ccataatctt taatacttcc taaactttct taataagaga agctctattc ctgacactac 120
ctctcatttg caagggtcaaa tcatcattag ttttgtagtc tattaactgg gtttgcttag 180
gtcaggcatt attattacta accttattgt taatattcta accataagaa ttaaactatt 240
aatggtgaat agagtttttc actttaacat aggcctatcc cactgggtggg atacgagcca 300
attcgaaaga aaagtcagtc atgtgctttt cagaggatga aagcttaaga taaagactaa 360
aagtgtttga tgctggaggt gggagtggtta ttatataggt ctgagccaag acatgtgata 420
atcactgtag tagtagctgg aaagagaaat ctgtgactcc aattagccag ttcctgcaga 480
ccttgtaggg actagaggaa gaatgctcct ggctgttttg tactgcctgc tgtggagttt 540
ccagacctcc gctggccatt tccctagagc ctgtgtctcc tctaagaacc tgatggagaa 600
ggaatgctgt ccaccgtgga gcggggacag gagtccctgt ggccagcttt caggcagagg 660
ttcctgtcag aatatacttc tgtccaatgc accacttggg cctcaatttc ccttcacagg 720
ggtggatgac cgggagtcgt ggccttccgt cttttataat aggacctgcc agtgctctgg 780
caacttcatt ggattcaact gtggaaactg caagtttggc ttttggggac caaactgcac 840
agagagacga ctcttggtga gaagaaacat cttcgatttg agtgccccag agaaggacaa 900
atTTTTTgcc tacctcactt tagcaaagca taccatcagc tcagactatg tcatcccat 960
agggacctat ggccaaatga aaaatggatc aacacccatg ttttaacgaca tcaatattta 1020
tgacctcttt gtctggatgc attattatgt gtcaatggat gcactgcttg ggggatctga 1080
aatctggaga gacattgatt ttgcccataga agcaccagct tttctgcctt ggcatagact 1140

cttcttgttg cggtgggaac aagaaatcca gaagctgaca ggagatgaaa acttcactat 1200
 tccatattgg gactggcggg atgcagaaaa gtgtgacatt tgcacagatg agtacatggg 1260
 aggtcagcac cccacaaatc ctaacttact cagcccagca tcattcttct cctcttggca 1320
 gattgtctgt agcogattgg aggagtacaa cagccatcag tctttatgca atggaacgcc 1380
 cgagggacct ttaaggcgta atcctggaaa ccatgacaaa tccagaacct caaggctccc 1440
 ctcttcagct gatgtagaat ttgcttgag ttgacccaa tatgaatctg gttccatgga 1500
 taaagctgcc aatttcagct ttagaaatac actggaagga ttgctagtc cacttactgg 1560
 gatagcggat gcctctcaaa gcagcatgca caatgccttg cacatctata tgaatggaac 1620
 aatgtcccag gtacagggat ctgccaacga tcctatcttc cttcttcacc atgcatttgt 1680
 tgacagtatt ttgagcagt ggctccgaag gcaccgtcct cttcaagaag tttatccaga 1740
 agccaatgca ccattggac ataaccggga atcctacatg gttcctttta taccactgta 1800
 cagaaatggg gatttcttta tttcatccaa agatctgggc tatgactata gctatctaca 1860
 agattcagac ccagactctt ttcaagacta cattaagtcc tatttggaac aagcgagtcg 1920
 gatctggtea tggctccttg gggcggcgat ggtaggggcc gtctcactg cctgctggc 1980
 agggcttggt agcttgctgt gtcgtcacia gagaaagcag cttcctgaag aaaagcagcc 2040
 actcctcatg gagaaagagg attaccacag cttgtatcag agccatttat aaaaggctta 2100
 ggcaatagag tagggccaaa aagcctgacc tcactctaac tcaaagtaat gtccagggtc 2160
 ccagagaata tctgctggta tttttctgta aagaccattt gcaaaattgt aacctaatc 2220
 aaagtgtagc cttcttccaa ctcaggtaga acacacctgt ctttgtcttg ctgttttcac 2280
 tcagcccttt taacattttc ccctaagccc atatgtctaa ggaaaggatg ctatttggt 2340
 atgaggaact gttatttgta tgtgaattaa agtgccttta tttt 2384

<210> 23

<211> 30

<212> DNA

<213> Artificial

<220>

<223> Artificial sequence is a Probe

<400> 23

cctgaggcac ctggaaggag gctgcagtgt

30

<210> 24

<211> 27
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Probe

<400> 24
atcaaaaaac aagcatggcc tcacacc

27

<210> 25
<211> 22
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 25
ctgcttcgct gcatcgctga aa

22

<210> 26
<211> 22
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 26
cagactcctc cagtcaggta ca

22

<210> 27
<211> 21
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 27
gcttggtggt taaaacttac c

21

<210> 28
<211> 20
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 28
tgaacagttc tgttggtgta

20

<210> 29
<211> 27
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Probe

<400> 29
ctgcctgcct atgtgacgac aatccgg

27

<210> 30
<211> 20
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 30
caattttggt ggagaacccg

20

<210> 31
<211> 20
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 31
gctgtcggag gtatatggtg

20

<210> 32
<211> 28
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Probe

<400> 32
catttcagag agtaacatgg actacaca

28

<210> 33
<211> 21
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 33
gcaagtgccca atgatcagag g

21

<210> 34
<211> 23
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Primer

<400> 34
atatagactc aggtatacac act

23

<210> 35
<211> 27
<212> DNA
<213> Artificial

<220>
<223> Artificial sequence is a Probe

<400> 35
tcccatcaga atccaaacaa gaggaag

27

CLAIMS:

1. A method of conducting an intra-operative molecular diagnostic assay comprising the steps of: analyzing a lymph node tissue sample from a patient by extracting RNA from the lymph node; amplifying and detecting the RNA; and determining if the presence of more than one Marker exceeds a cut-off value, wherein said RNA is extracted in less than 8 minutes and wherein the method is conducted within a period of no more than 35 minutes.
2. The method of claim 1 wherein the said RNA is extracted in less than 6 minutes.
3. The method of claim 1 for use in detecting micrometastasis.
4. The method of claim 1 for use in detecting breast cancer metastasis.
5. The method of any one of claims 1 to 3 wherein the Markers detect the expression of a gene corresponding to Seq ID 17 and to one or more of Seq ID Nos 14-16 and Seq ID No 18.
6. The method of any one of claims 1 to 5 wherein nucleic acid amplification and detection is conducted by PCR.
7. The method of claim 6 wherein said PCR is RT-PCR.
8. The method of claim 7 wherein the reverse transcription reaction is conducted for less than 9 minutes.
9. The method of claim 8 wherein the reverse transcription reaction is conducted for about 3 minutes.
10. The method of any one of claims 7 to 9, wherein one or more internal control reagents is added to the reverse transcription reaction.
11. The method of any one of claims 1 to 10 wherein RNA is extracted from the

lymph node tissue sample by:

- (a) lysing the cells containing the RNA and producing a homogenate of them,
- (b) optionally, diluting the homogenate,
- (c) contacting the wetted, homogenized working sample with a substrate containing, or to which is affixed, a material to which RNA binds,
- (d) allowing the sample to bind to the substrate,
- (e) removing contaminants and interferents,
- (f) drying the substrate, and
- (g) eluting RNA from the substrate; and

employing filtration in extraction steps.

12. A method of conducting an intra-operative molecular diagnostic assay for use in detecting melanoma metastasis comprising the steps of: analyzing a lymph node tissue sample from a patient by extracting RNA from the lymph node; amplifying and detecting the RNA; and determining if the presence of more than one Marker exceeds a cut-off value, wherein said RNA is extracted in less than 8 minutes and wherein the method is conducted within a period of no more than 35 minutes.

13. The method of claim 12 wherein the said RNA is extracted in less than 6 minutes.

14. The method of claim 12 for use in detecting micrometastasis.

15. The method of claim 14 wherein the Markers are Seq ID No 21 and Seq ID No 22.

16. The method of any one of claims 1 to 4 wherein nucleic acid amplification and detection is conducted by PCR.


17. The method of claim 16 wherein said PCR is RT-PCR.

18. The method of claim 17 wherein the reverse transcription reaction is conducted for less than 9 minutes.

19. The method of claim 18 wherein the reverse transcription reaction is conducted for about 3 minutes.
20. The method of any one of claims 17 to 19, wherein one or more internal control reagents is added to the reverse transcription reaction.
21. The method of any one of claims 12 to 20 wherein RNA is extracted from the lymph node tissue sample by:
 - (a) lysing the cells containing the RNA and producing a homogenate of them,
 - (b) optionally, diluting the homogenate,
 - (c) contacting the wetted, homogenized working sample with a substrate containing, or to which is affixed, a material to which RNA binds,
 - (d) allowing the sample to bind to the substrate,
 - (e) removing contaminants and interferents,
 - (f) drying the substrate, and
 - (g) eluting RNA from the substrate; andemploying filtration in extraction steps.


1/5

Figure 1




2/5

Figure 2




Expression of Mammaglobin in Different Breast Cell Types



4/5

Figure 4



5/5

Figure 5

