



(51) International Patent Classification:
C12M 1/00 (2006.01) C12M 1/36 (2006.01)

(21) International Application Number:
PCT/EP2020/053654

(22) International Filing Date:
12 February 2020 (12.02.2020)

(25) Filing Language: English

(26) Publication Language: English

(30) Priority Data:
19158305.3 20 February 2019 (20.02.2019) EP

(71) Applicant: SARTORIUS STEDIM BIOTECH GMBH
[DE/DE]; August-Spindler-Strasse 11, 37079 Goettingen
(DE).


(72) Inventors: MATUSZCZYK, Jens-Christoph; c/o Sartorius Stedim Biotech GmbH, August-Spindler-Strasse 11, 37079 Goettingen (DE). SCHOLZ, Jochen; c/o Sartorius Stedim Biotech GmbH, August-Spindler-Strasse 11, 37079 Goettingen (DE). ZIJLSTRA, Gerben; c/o Sartorius Stedim Biotech GmbH, August-Spindler-Strasse 11, 37079 Goettingen (DE). GRELLER, Gerhard; c/o Sartorius Stedim Biotech GmbH, August-Spindler-Strasse 11, 37079 Goettingen (DE).

(74) Agent: PRINZ & PARTNER MBB PATENT- UND RECHTSANWÄLTE; Rundfunkplatz 2, 80335 München (DE).

(81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DJ, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN,

(54) Title: BIOPROCESS DEVICE ASSEMBLY AND INOCULATION METHOD

Fig. 1



(57) Abstract: A bioprocess device (10) assembly comprises at least a first bioreactor (12), a second bioreactor (14), a transfer system (18) for transferring a certain volume from the first bioreactor (12) into the second bioreactor (14). A biomass sensor (16) is associated to the first bioreactor (12). A control unit (24) is configured for receiving measurement data from the biomass sensor (16) and for controlling the transfer system (18). An inoculation method makes use of such a device assembly (10) by following steps: the biomass sensor (16) measures at least one parameter of the biomass contained in the first bioreactor (12); the biomass sensor (16) sends measurement data to the control unit (24); and the control unit (24) evaluates the data received from the biomass sensor (16) and decides whether the viable cell density is sufficient for inoculation, determines an inoculation volume, and controls the transfer from the first bioreactor (12) into the second bioreactor (14) by transfer system (18).

HR, HU, ID, IL, IN, IR, IS, JO, JP, KE, KG, KH, KN, KP, KR, KW, KZ, LA, LC, LK, LR, LS, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SA, SC, SD, SE, SG, SK, SL, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, WS, ZA, ZM, ZW.

(84) Designated States (*unless otherwise indicated, for every kind of regional protection available*): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, ST, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, KM, ML, MR, NE, SN, TD, TG).

Published:

— with international search report (Art. 21(3))

Bioprocess Device Assembly and Inoculation Method

The invention relates to a bioprocess device assembly. The invention further relates to an inoculation method performed within a bioprocess.

5 Process intensification usually refers to a strategy aimed at transforming established processes into more efficient and economical processes. The overall productivity is increased by streamlining operations and applying more advanced technologies. In the field of industrial bioprocessing it is known to intensify fed-batch processes for mammalian cells by means of high cell seeding by an approach relying on a seed bioreactor (N-1 stage) operated in perfusion mode in
10 order to reach elevated cell densities. Subsequently, the fed-batch bioreactor (N stage) is started at a cell concentration which is substantially higher than what is done with conventional fed-batch processes (Stettler et al: "Intensification of large-scale mammalian fed-batch processes", Pharmaceutical Bioprocessing, 2018 Volume 6, Issue 1).

15 Nevertheless, transfer of the inoculum from one bioreactor to the other is complicated and involves certain risks. In order to determine the viable cell density (VCD), which is one of the main criteria for a "good" inoculum" and used here as the representative method for measuring the viable biomass in a bioreactor, a sample has to be taken manually and examined off-line. If the measured VCD is
20 not sufficient for effective inoculation, another sample has to be taken and examined later, and so on. If the intervals between the sample takings are long, e.g. one day, process time is usually wasted as it is very unlikely that a sample is taken exactly at the time when the desired, pre-determined VCD threshold value is first reached. On the other hand, if samples are taken frequently, this means
25 increased effort for the process operator.

Moreover, after calculating a necessary inoculation volume depending on the actual VCD, transfer of the correct volume is controlled manually, which is also a laborious and error-prone task.

In general, when samples are taken and medium is transferred, each physical interaction with the bioreactors and the medium increases the risk of contamination.

It is an object of the invention to offer more control over a bioprocess, especially
5 with regard to a precise transfer of cell cultures for inoculation, and to make the bioprocess thus more efficient.

The above object is achieved by a bioprocess device assembly according to claim 1 and also by an inoculation method performed within a bioprocess according to claim 6. Advantageous and expedient embodiments of the invention are
10 apparent from the associated dependent claims.

The invention provides a bioprocess device assembly, comprising at least a first bioreactor and a second bioreactor, and a transfer system for transferring a certain volume of a medium contained in the first bioreactor into the second bioreactor. The bioprocess device assembly further comprises a biomass sensor
15 associated to the first bioreactor for measuring at least one parameter relating to the quantity and/or quality of the biomass contained in the first bioreactor, and a control unit configured for receiving measurement data from the biomass sensor and controlling the transfer system. The control unit is parametrized such that, based on an evaluation of the measurement data received from the biomass
20 sensor, it decides whether the biomass status, preferably the viable cell density (VCD), is sufficient for inoculation, determines a necessary inoculation volume, and controls the transfer system to transfer the determined inoculation volume from the first bioreactor into the second bioreactor.

The invention is based on the general finding that online-monitoring of critical
25 bioprocess parameters, especially VCD, offers the potential of significantly improving process stability, reducing process time and enhancing process performance. In fact, VCD of the inoculum is considered as decisive parameter for quality and reproducibility of subsequent steps. Accordingly, the most important requirement of optimal control of all follow up steps is inoculation with a precise
30 VCD.

According to the invention, a biomass sensor is provided for online measuring of at least one parameter relating to the quantity and/or quality of the biomass

contained in the first (pre-culture) bioreactor. Online monitoring of this parameter, especially VCD, makes it possible to perform inoculation as soon as the parameter has reached a desired, pre-determined threshold value, i.e. without substantial delay.

5 The invention is based on the further finding that VCD can be determined online by measuring the capacitance of the biomass when the corresponding bioreactor is operated in a perfusion mode. In this case, certain cell parameters that would otherwise change significantly over time, e.g. increasing cell diameter due to swelling of the cells, remain rather constant so that their influence on the
10 capacitance-VCD correlation is not detrimental.

As already indicated, it is preferred that the biomass sensor is at least capable of measuring capacitance. This can be either capacitance measurements at a single frequency as well as multi-frequency capacitance measurements (also referred to as capacitance scanning or capacitance spectroscopy).

15 It is to be noted that the quantity and/or quality of the biomass can be measured by a "soft" biomass sensor, i. e. by measuring a parameter or a combination of parameters which is/are dependent on the biomass. Therefore, in addition or in the alternative, the biomass sensor can also be capable of measuring at least one of the following: electrical impedance; oxygen transmission rate (OTR); oxygen
20 uptake rate (OUR); scattered light; radiation intensity; concentration of a metabolite (e.g. glucose); metabolite (e.g. glucose) consumption rate; product formation rate. Regarding the latter, spectrometry allows monitoring the protein content in the harvest stream. Together with a constant cell-specific product formation rate this can be used as a reliable indicator for biomass. If necessary or expedient, several
25 parameters can be measured and evaluated in combination.

The biomass sensor and the control unit are preferably configured for monitoring the at least one parameter and controlling the transfer system in real-time. This means that measurement data is delivered and evaluated immediately after collection, so that the transfer system can be activated without (significant)
30 delay after certain conditions are met.

The transfer system for transferring medium from the first bioreactor into the second bioreactor can be of any suitable type. In particular, the transfer system

can include at least one of the following controllable conveyor devices: pump; pneumatic drive; gravimetric feeder. A flow meter and/or other devices can be used to control the transferred inoculation volume.

5 The invention further provides an inoculation method performed within a bioprocess making use of a device assembly according to the invention. The inoculation method comprises the following steps: the biomass sensor constantly or repeatedly measures at least one parameter relating to the quantity and/or quality of the biomass contained in the first bioreactor; the biomass sensor sends measurement data to the control unit; and the control unit evaluates the
10 measurement data received from the biomass sensor and, based on an evaluation of the measurement data decides whether the biomass status, preferably the viable cell density (VCD), is sufficient for inoculation, determines a necessary inoculation volume, and controls the transfer system to transfer the determined inoculation volume from the first bioreactor into the second bioreactor.

15 Preferably, the steps of evaluating the measurement data received from the biomass sensor and controlling the transfer system to transfer the determined inoculation volume from the first bioreactor into the second bioreactor are performed in real-time.

To ensure optimum process performance, the necessary inoculation volume
20 should be calculated based on the evaluated measurement data of the biomass sensor, because the required inoculation volume largely depends on the viable cell density of the inoculum.

Transfer of the inoculation volume can basically be performed in many different ways. In particular, one of the following driving principles can be used: hydraulic
25 pressure; pneumatic pressure; gravity.

The first bioreactor is preferably operated in a perfusion mode. Perfusion is not only beneficial in view of bioprocess intensification, but in perfusion mode the correlation between the measured parameter, especially capacitance, is not impaired by side effects like swelling of the cells etc.

30 Apart from the calculation of the necessary inoculation volume, its advantageous to adapt or adjust the initial conditions in the second bioreactor and/or process control and process parameters of the second bioreactor

depending on at least one of the following: amount or concentration of a metabolite in the transferred inoculation volume; one or more physicochemical characteristics of the transferred inoculation volume, in particular pH, oxygen saturation, temperature.

- 5 The inventive concept is scalable. If a cascade of pre-culture bioreactors (N₁, ... N-1) is used in the bioprocess, the method according to the invention can be performed multiple times within the cascade of bioreactors.

Further features and advantages of the invention will become apparent from the following description and from the accompanying drawings to which reference
10 is made. In the drawings:

- Figure 1 shows an example setup of a bioprocess device assembly according to the invention using a control system;
- Figure 2 shows an example correlation between viable cell density and capacitance; and
- 15 - Figure 3 shows generally a cascade of a plurality of controlled bioreactors in an assembly according to the invention.

Figure 1 schematically depicts a possible setup of a bioprocess device assembly 10 for performing an inoculation within a fed-batch bioprocess. The bioprocess device assembly 10 includes at least two bioreactors 12, 14. The first
20 bioreactor 12 is used as a seed bioreactor (N-1 stage) and preferably operated in a perfusion mode, while a second bioreactor 14 is used as a fed-batch bioreactor (N stage). The bioreactors 12, 14 shown in Figure 1 are both rocking motion (RM) bioreactors.

The first bioreactor 12 is equipped with a biomass sensor 16 for measuring at
25 least one parameter relating to the quantity and/or quality of the biomass contained in the first bioreactor 12. The biomass sensor 16 preferably is or includes a capacitance sensor. In addition, or alternatively, the biomass sensor 16 may be capable of measuring at least one of the following parameters or any derivative thereof: electrical impedance; oxygen transmission rate (OTR); oxygen uptake rate
30 (OUR); scattered light; radiation intensity; concentration of a metabolite (e.g. glucose); metabolite (e.g. glucose) consumption rate; product formation rate.

A transfer system 18 makes it possible to transfer a medium contained in the first bioreactor 12 into the second bioreactor 14. In the example embodiment shown in Figure 1 the transfer system 18 includes a hose line 20 connecting the two bioreactors 12, 14 and a controllable pump 22 integrated in the hose line 20.

5 However, other automatically or semi-automatically controllable transfer concepts may be employed. In particular, the transfer system 18 may also include a pneumatic drive, a gravimetric feeder, or the like.

The biomass sensor 16 and the transfer system 18 are connected to a control unit 24 which is configured for receiving measurement data from the biomass
10 sensor 16 and for controlling the transfer system 18. The control unit 24 can be a specific unit mainly dedicated for these purposes, or it can be a flexible digital control unit which is generally capable of performing a variety of process control operations. In any event, the control unit 24 may be part of, or connected to, a higher level control system for controlling several or all steps in a bioprocess.

15 The control unit 24 is capable of controlling the transfer system 18 such that a previously determined amount of medium contained in the first bioreactor 12, in particular an inoculation volume, is transferred from the first bioreactor 12 into the second bioreactor 14. For exactly measuring the weight or volume of the transferred amount of medium one or more weight and/or volume sensors
20 associated to one or both of the bioreactors 12, 14 are connected to the control unit 24.

The control unit 24 itself is parametrized such that, based on an evaluation of the measurement data received from the biomass sensor 16 in the first bioreactor 12, the control unit 24 performs the following steps:

- 25 (a) deciding whether the viable cell density (VCD) is sufficient for inoculation,
(b) determining a necessary inoculation volume, and
(c) controlling the transfer system 18 to transfer the determined inoculation volume from the first bioreactor 12 into the second bioreactor 14.

The parametrization of the control unit 24 is based on a known correlation
30 between VCD and the parameter(s) measured by the biomass sensor 16, especially capacitance. The correlation was determined earlier empirically or,

especially if empirical data is not available, in parallel to the actual seed steps. The correlation can also be obtained by modelling. The correlation is stored in the control unit 24. An example of a substantially linear correlation between VCD and capacitance is shown in the diagram of Figure 2.

5 The inoculation within the bioprocess is performed as follows. The biomass sensor 16 continuously or repeatedly provides measurement data to the control unit 24, and the control unit 24 continuously or repeatedly evaluates the measurement data in real-time (online monitoring). When the control unit 24 determines, based on the above-mentioned known correlation, that the VCD is
10 sufficient or has reached a desired, pre-determined threshold value, inoculation is immediately and automatically initiated.

Before the actual medium transfer from the first bioreactor 12 into the second bioreactor 14, the control unit 24 determines a necessary amount of biomass for inoculation (inoculation volume) according to the present circumstances, especially
15 based on the measured VCD. The medium transfer is then automatically activated by the control unit 24.

In the example embodiment shown in Figure 1 the transfer of the inoculation volume is mainly performed by operating the pump 22 in the hose line 20 between the two bioreactors 12, 14. Operation of the pump 22 (at least starting and stopping
20 the pump 22) is controlled by the control unit 24. Additionally, valves for opening and closing the bioreactor ports connected to the hose line 20 may also be controlled by the control unit 24. If another type of transfer system is used, e.g. based on a pneumatic drive, a gravimetric feeder, or the like, the control unit 24 controls these components accordingly.

25 As soon as the calculated amount of inoculation volume is transferred by the transfer system 18, the control unit 24 automatically terminates the medium transfer based on the measurement data received from the weight and/or volume sensor(s).

30 The amount or concentration of metabolites, e.g. glucose, contained in the inoculation volume and/or one or more physicochemical characteristics of the inoculation volume, such as pH, oxygen saturation (dissolved oxygen, DO), temperature etc., are determined by calculation or measurement. Based on the

determined value(s), the initial conditions in the second bioreactor 14 are adapted or adjusted accordingly to the new amount of metabolites and/or physicochemical characteristics. The same can be done with the process control and process parameters (PID) of the second bioreactor 14, i.e. they can be adapted or adjusted
5 according to the determined value(s).

An optional biomass sensor may be provided in the second bioreactor 14 for examining the effectively transferred biomass.

The above-described concept, or parts of it, can be scaled up as indicated by the dashed lines in Figure 1. In particular, the transfer of medium between several
10 pre-culture bioreactors (N-x stage to N-x+1 stage) can be automatically initiated and automatically performed as described above with the same control unit 24 or several control units. A corresponding cascade of controlled bioreactors is generally shown in Figure 3.

15

List of Reference Signs

	10	bioprocess device assembly
	12	first bioreactor
	14	second bioreactor
5	16	biomass sensor
	18	transfer system
	20	hose line
	22	pump
	24	control unit
10		

Claims

1. A bioprocess device assembly (10), comprising
at least a first bioreactor (12) and a second bioreactor (14),
a transfer system (18) for transferring a certain volume of a medium contained
5 in the first bioreactor (12) into the second bioreactor (14),
a biomass sensor (16) associated to the first bioreactor (12) for measuring at
least one parameter relating to the quantity and/or quality of the biomass contained
in the first bioreactor (12), and
a control unit (24) configured for receiving measurement data from the biomass
10 sensor (16) and controlling the transfer system (18),
the control unit (24) being parametrized such that, based on an evaluation of
the measurement data received from the biomass sensor (16), the control unit
(a) decides whether the biomass status, preferably the viable cell density, is
sufficient for inoculation,
15 (b) determines a necessary inoculation volume, and
(c) controls the transfer system (18) to transfer the determined inoculation
volume from the first bioreactor (12) into the second bioreactor (14).
2. The bioprocess device assembly (10) according to claim 1, characterized
in that the biomass sensor (16) is capable of measuring capacitance.
- 20 3. The bioprocess device assembly (10) according to claim 1 or 2,
characterized in that the biomass sensor (16) is capable of measuring at least one
of the following: electrical impedance; oxygen transmission rate; oxygen uptake
rate; scattered light; radiation intensity; concentration of a metabolite; metabolite
consumption rate; product formation rate.
- 25 4. The bioprocess device assembly (10) according to any of the preceding
claims, characterized in that the biomass sensor (16) and the control unit (24) are
configured for monitoring the at least one parameter and controlling the transfer
system (18) in real-time.
- 30 5. The bioprocess device assembly (10) according to any of the preceding
claims, characterized in that the transfer system (18) includes at least one of the

following controllable conveyor devices: pump (22); pneumatic drive; gravimetric feeder.

6. An inoculation method performed within a bioprocess, the method making use of a device assembly (10) according to any of the preceding claims and comprising the following steps:

- the biomass sensor (16) constantly or repeatedly measures at least one parameter relating to the quantity and/or quality of the biomass contained in the first bioreactor (12);

- the biomass sensor (16) sends measurement data to the control unit (24); and

- the control unit (24) evaluates the measurement data received from the biomass sensor (16) and, based on an evaluation of the measurement data

- (a) decides whether the biomass status, preferably the viable cell density, is sufficient for inoculation,

- (b) determines a necessary inoculation volume, and

- (c) controls the transfer system (18) to transfer the determined inoculation volume from the first bioreactor (12) into the second bioreactor (14).

7. The method according to claim 6, characterized in that the steps of evaluating the measurement data received from the biomass sensor (16) and controlling the transfer system (18) to transfer the determined inoculation volume from the first bioreactor (12) into the second bioreactor (14) are performed in real-time.

8. The method according to claim 6 or 7, characterized in that the necessary inoculation volume is calculated based on the evaluated measurement data of the biomass sensor (16).


9. The method according to any of the preceding claims, characterized in that the transfer of the inoculation volume is performed based on one of the following driving principles: hydraulic pressure; pneumatic pressure; gravity.

10. The method according to any of claims 6 to 9, characterized in that the first bioreactor (12) is operated in a perfusion mode.



11. The method according to any of claims 6 to 10, characterized in that the initial conditions in the second bioreactor (14) and/or process control and process parameters of the second bioreactor (14) are adapted or adjusted depending on at least one of the following: amount or concentration of a metabolite in the transferred inoculation volume; one or more physicochemical characteristics of the transferred inoculation volume, in particular pH, oxygen saturation, temperature.

12. The method according to any of claims 6 to 11, characterized in that the method is performed multiple times with a cascade of pre-culture bioreactors.

Fig. 1



2 / 2

Fig. 2**Fig. 3**

INTERNATIONAL SEARCH REPORT

International application No
PCT/EP2020/053654

A. CLASSIFICATION OF SUBJECT MATTER
INV. C12M1/00 C12M1/36
ADD.

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)
C12M

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

EPO-Internal, WPI Data

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	WO 2015/003012 A2 (BERTEAU OLIVIER [US]; YOON SEONGKYU [US] ET AL.) 8 January 2015 (2015-01-08) page 3, line 2 - line 24 page 12, line 2 - line 12 page 15, line 12 - line 25 page 17, line 17 - line 31 page 20, line 1 - line 10 page 20, line 31 - page 21, line 21 page 21, line 23 - line 30 page 22, line 5 - line 10 page 23, line 20 - line 24 page 30, line 31 - page 31, line 3 page 44, line 20 - page 45, line 7 figures 2-9e ----- -/-	1-12



Further documents are listed in the continuation of Box C.



See patent family annex.

* Special categories of cited documents :

"A" document defining the general state of the art which is not considered to be of particular relevance

"E" earlier application or patent but published on or after the international filing date

"L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)

"O" document referring to an oral disclosure, use, exhibition or other means

"P" document published prior to the international filing date but later than the priority date claimed

"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention

"X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone

"Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art

"&" document member of the same patent family

Date of the actual completion of the international search

5 March 2020

Date of mailing of the international search report

13/03/2020

Name and mailing address of the ISA/

European Patent Office, P.B. 5818 Patentlaan 2
NL - 2280 HV Rijswijk
Tel. (+31-70) 340-2040,
Fax: (+31-70) 340-3016

Authorized officer

Panzica, Gian Paolo

INTERNATIONAL SEARCH REPORT

International application No

PCT/EP2020/053654

C(Continuation). DOCUMENTS CONSIDERED TO BE RELEVANT		
Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	WO 2018/115161 A1 (HOFFMANN LA ROCHE [CH]; HOFFMANN LA ROCHE [US]) 28 June 2018 (2018-06-28) page 7, line 14 - line 24 page 8, line 16 - line 20 page 8, line 26 - line 30 page 9, line 8 - line 20 page 26, line 21 - line 27 page 34, line 28 - page 35, line 10 figures 1, 3, 4 -----	1-12
A	WO 2006/071716 A2 (BIOGEN IDEC INC [US]; WEBB MARCUS [US]) 6 July 2006 (2006-07-06) page 2, line 7 - line 9 page 4, line 17 - line 27 page 7, line 26 - line 27 page 9, line 17 - line 20 -----	2

INTERNATIONAL SEARCH REPORT

Information on patent family members

International application No

PCT/EP2020/053654

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 2015003012 A2	08-01-2015	US 2016145563 A1 WO 2015003012 A2	26-05-2016 08-01-2015
WO 2018115161 A1	28-06-2018	CN 110023482 A EP 3559203 A1 JP 2020504601 A US 2020063086 A1 WO 2018115161 A1	16-07-2019 30-10-2019 13-02-2020 27-02-2020 28-06-2018
WO 2006071716 A2	06-07-2006	AU 2005322159 A1 CA 2593374 A1 CN 101443444 A EP 1844138 A2 JP 2008526203 A US 2009104594 A1 WO 2006071716 A2	06-07-2006 06-07-2006 27-05-2009 17-10-2007 24-07-2008 23-04-2009 06-07-2006