



- (51) **International Patent Classification:**  
*A61B 17/42* (2006.01) *A61B 18/04* (2006.01)
- (21) **International Application Number:**  
PCT/US2015/032009
- (22) **International Filing Date:**  
21 May 2015 (21.05.2015)
- (25) **Filing Language:** English
- (26) **Publication Language:** English
- (30) **Priority Data:**  
62/002,082 22 May 2014 (22.05.2014) US
- (71) **Applicant:** AEGEA MEDICAL INC. [US/US]; 2686 Middlefield Road, Suite A, Redwood City, CA 94063 (US).
- (72) **Inventors:** CHEE, Uriel, Hiram; 2686 Middlefield Road, Suite A, Redwood City, CA 94063 (US). PELIKS, Robert, Bilgor; 2686 Middlefield Road, Suite A, Redwood City, CA 94063 (US). MAGEN, Hugh, Edward; 2686 Middlefield Road, Suite A, Redwood City, CA 94063 (US). GURSKIS, Donnell, William; 2686 Middlefield Road, Suite A, Redwood City, CA 94063 (US). BACICH, Steven, Robert; 2686 Middlefield Road, Suite A, Redwood City, CA 94063 (US).
- (74) **Agents:** THOMAS, Justin et al.; Shay Glenn LLP, 2755 Campus Drive, Suite 210, San Mateo, CA 94403 (US).
- (81) **Designated States** (*unless otherwise indicated, for every kind of national protection available*): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IR, IS, JP, KE, KG, KN, KP, KR, KZ, LA, LC, LK, LR, LS, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SA, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.
- (84) **Designated States** (*unless otherwise indicated, for every kind of regional protection available*): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, ST, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, KM, ML, MR, NE, SN, TD, TG).
- Published:**  
— with international search report (Art. 21(3))

(54) **Title:** SYSTEMS AND METHODS FOR PERFORMING ENDOMETRIAL ABLATION




Fig. 1D

(57) **Abstract:** A method and system of providing therapy to a patient's uterus is provided, which can include any number of features. The method can include the steps of inserting a uterine device into the uterus and performing a uterine integrity test to determine that the uterus is intact and not perforated. If it is determined that the uterus is not perforated, a patency test can be performed to determine that the uterine device is not clogged or embedded in tissue. If the uterus is intact and the device is not clogged or embedded in tissue, the uterus can be treated with the uterine device, e.g., uterine ablation. Systems for performing these methods are also disclosed.



## **SYSTEMS AND METHODS FOR PERFORMING ENDOMETRIAL ABLATION**

### **CROSS REFERENCE TO RELATED APPLICATIONS**

- 5   **[0001]**    This application claims the benefit of US Provisional Application No. 62/002,082, filed May 22, 2014, titled "Systems and Methods for Performing Endometrial Ablation", which is incorporated herein by reference.

### **INCORPORATION BY REFERENCE**

- 10   **[0002]**    All publications and patent applications mentioned in this specification are herein incorporated by reference in their entirety to the same extent as if each individual publication or patent application was specifically and individually indicated to be incorporated by reference.

### **FIELD**

- 15   **[0003]**    The present disclosure generally relates to uterine procedures incorporating a distension media such as a fluid or a gas that could be used with endoscopic procedures or other visualization systems such ultrasound or fluoroscopy. The present disclosure is particular suited for endometrial ablation of the uterine lining. More specifically, the present disclosure relates to endometrial ablation with a heated vapor.

20

### **BACKGROUND**

- [0004]**    Endometrial ablation (i.e., the removal or destruction of the endometrial lining of the uterus) is used as an alternative to hysterectomy for treating menorrhagia, or other uterine diseases. One prior technique for performing endometrial ablation employs a resectoscope (i.e.,  
25   a hysteroscope with a built-in wire loop or other ablative devices) that is inserted transcervically into the uterus, and uses radio-frequency electrical current (RF current) to remove or coagulate the endometrial tissue. These standard techniques typically are performed in a hospital setting and importantly utilize hysteroscopy for visualization of the procedure while treating the uterine lining.

- 30   **[0005]**    Some approaches make use of heated fluid to ablate the endometrium. For example, early journal articles describe the use of steam to treat uterine hemorrhage. The use of steam for this purpose was later discredited, apparently due to patient morbidity and mortality. See, e.g., Fuller U.S. Patent No. 6,139,571. More recent descriptions of the use of injecting hot fluid into the uterus have been described. Uterine therapies employing a contained fluid have also been  
35   described.

[0006] In an effort to simplify the procedure, approaches have been developed that do not require concurrent hysteroscopic visualization. In practice, many of these techniques recommend that the physician or user employ hysteroscopy to visualize and inspect the uterine cavity prior to performing the endometrial ablation procedure. In addition, hysteroscopy may be employed at the conclusion of the endometrial ablation procedure as a method to inspect the uterine cavity post treatment. During this hysteroscopic inspection, the physician is verifying that the uterine cavity is not perforated although perforations may not be readily apparent even with hysteroscopic visualization. In general, a physician seeks to avoid perforations for many reasons including the potential for unintended injuries to neighboring organs and maintaining or confining the treatment area to specifically the uterine cavity in the case of endometrial ablation procedures.

[0007] Endometrial ablation techniques that do not require active hysteroscopic visualization during treatment operation are commonly referred to as “blind” techniques since the physician is using tactile feel, or markers and indicia on the endometrial ablation device to indicate proper placement of the device in the uterine cavity. One of these particular devices utilizes a balloon-based system using heated saline as the thermal energy source for the ablation of tissue. High frequency, or radiofrequency (RF), energy has also been used to perform thermal ablation of endometrial tissue. Current products for performing endometrial ablation include the NOVASURE<sup>®</sup> procedure and a system marketed under the trade name THERMACHOICE<sup>®</sup>, by Ethicon, Inc. of Somerville, N.J. Cryogenic ablation, or “cryoablation,” such as HER OPTION<sup>®</sup> from American Medical Systems, Inc., is another endometrial treatment approach. All of the products above are characterized as “blind” or not requiring direct hysteroscopic visualization during the treatment.

[0008] In utilizing an endometrial ablation technology that does not require hysteroscopic visualization, it would be beneficial to employ a test to verify that the uterine cavity is intact or unperforated prior to performing the treatment. Such tests are referred to as uterine integrity tests and these tests can be performed with endometrial ablation procedures and any procedure of the uterus or hollow body cavity or organ. In addition, these tests can be used with hysteroscopic procedures since a perforation may not be readily detected even under direct vision.

[0009] Integrity tests employ saline or gas, preferably carbon dioxide gas, as agents to verify if the uterine cavity is intact in regards to holding fluid or gas pressure. The gas or fluid is supplied under pressure to the uterine cavity and a leak in the uterine cavity, whether it is a perforation, an unsealed cervical canal, or the effect of excess fluid exiting the fallopian tubes, can be discerned. Stern et al. (US 5,562,720) and Sampson et al. (US 6,554,780, US 6,743,184, US 6,872,183, and US 7,063,670) describe such pressure techniques while other approaches

check for fluid imbalances between an input source and output collection using volume measurements. Other approaches mention using flow rate and pressure measurements.

[0010] For monitoring the therapeutic energy application during the procedure, some technologies monitor the internal pressure of heated saline within a balloon that is placed within the uterus, or the impedance of radio frequency energy within the wall of the uterus. These technologies have automatic termination steps if the pressure, impedance, or volume levels reach certain thresholds. At the initiation of the procedures for all of the above mentioned systems, the inaccurate placement and management of the therapeutic device by the physician within the uterine cavity can reduce the ability to perform a safe and consistent ablation procedure. For these technologies, the ability of the device to perform a complete endometrial ablation procedure depends upon the tactile movements and placement of the device by the physician in terms of depth of placement or achieving complete contact with the interior lining with the delivery device mechanism. As an example, a radio frequency device deploys an RF array within the uterine cavity. The depth of insertion can vary by patient and physician. By not contacting the endometrium in more distal or proximal locations, the efficacy of the procedure could alter or be affected. For thermal balloon procedures, a similar effect can occur due to the depth of device placement. For both of the technologies, the system cannot provide information or assessment for the proper and consistent placement of the device in the uterine cavity. Alternatively, an amorphous ablation technology such as the free flow of heated saline does not have this physical limitation however even this technology, known as Hydro Thermo Ablation or HTA, relies heavily on the placement of the device by the physician to prevent leaking out of the endocervical canal.

[0011] Systems that depend heavily on physician manipulation or intervention during the procedure have been characterized as being "technique sensitive". The requirement for a physician to properly manipulate a delivery device and react appropriately before and during an endometrial ablation procedure can lead to an increase of adverse events or unreliable patient outcomes. This is particularly evident when threshold values prior to and during a procedure fall within a range requiring physician intervention or adjustment.

[0012] The following describe a control system that overcomes these deficiencies for a technology that provides vapor energy in combination with a number of unique and sequenced steps designed to ensure patient safety and procedural consistency. The control system also assesses delivery device placement and sealing within the uterine cavity and endocervical canal by providing information for the physician when additional pressure within a sealing balloon at the cervix may be needed.

## SUMMARY

[0013] An endometrial vapor ablation system, comprising a vapor generator, a uterine ablation device coupled to the vapor generator, the uterine ablation device configured for insertion into a uterine cavity of a patient, a control system coupled to the uterine ablation device and the vapor generator, the control system being configured to automatically perform a pre-procedural sterilization of the vapor generator and the uterine ablation device, the control system configured to prep the uterine ablation device for vapor delivery, the control system configured to test safety components of the uterine ablation device prior to and after insertion into the patient, the control system configured to monitor the application of vapor energy, monitor temperature reading at identified locations on the uterine ablation device, and the control system being configured to implement shutdown of the vapor generator following vapor delivery.

[0014] In some embodiments, the electronic controller is configured to prompt an end user of the uterine ablation device if the pre-procedural sterilization fails to initiate or complete.

[0015] In another embodiment, the electronic controller is configured to prompt an end user of the uterine ablation device if the testing of the safety components indicates an error in the system.

[0016] In some embodiments, the system further comprises a distal anchor balloon, a central sealing balloon, and a proximal positioning balloon disposed on a shaft of the uterine ablation device.

[0017] In one embodiment, the system further comprises a temperature sensor positioned near the central sealing balloon and configured to measure a temperature inside an endocervix of the patient.

[0018] In yet another embodiment, the system further comprises one or more pressure sensors configured to measure a pressure of each of the distal anchor balloon, the central sealing balloon, and the proximal positioning balloon.

[0019] A method of providing vapor therapy to a uterus of a patient is provided, comprising the steps of sensing a fluid level of a vapor generator, heating the vapor generator to prepare a condensable vapor, priming the condensable vapor through a uterine ablation device coupled to the vapor generator to sterilize the uterine ablation device, comparing a pressure sensor value within the uterine ablation device to a pressure sensor value within the vapor generator, verifying the integrity of a plurality of positioning and sealing balloons of the uterine ablation device, delivering condensable vapor from the uterine ablation device into the uterus, and during the delivering step, monitoring a temperature in a cervix of the patient and automatically terminating delivery of condensable vapor if the temperature exceeds a threshold value.

[0020] In some embodiments, the threshold value comprises 44 degrees C.

### BRIEF DESCRIPTION OF THE DRAWINGS

- [0021] Figs. 1A-1D illustrate one embodiment of a uterine ablation device.
- [0022] Fig. 2 shows an embodiment of a uterine ablation device inserted into a uterus.
- 5 [0023] Fig. 3 illustrates an integrity test of the uterine ablation device.
- [0024] Fig. 4 illustrates a patency test of the uterine ablation device.
- [0025] Fig. 5 charts the relationship between outflow valve pressure, uterine cavity pressure, and flow through a uterine ablation device.

### DETAILED DESCRIPTION

10 [0026] Fig. 1A illustrates a uterine ablation device 100 sized and configured to access the endometrium of a uterus and to deliver a heated vapor to the uterus to ablate uterine tissue. The device can be configured to ablate and treat the endometrial lining of the uterus as an alternative to hysterectomy for treating menorrhagia or other uterine diseases. In some embodiments, the

15 device 100 can be configured to gain access to the uterus by being inserted through a cannula or hysteroscope. The device 100 can include shaft 102, handle 104, distal tip 106, vapor ports 107, distal anchor or distal balloon 108, central or sealing balloon 110, proximal or positioning balloon 112, and connection lumens 118, which can couple the uterine ablation device to a control system (not shown) comprising a computer, a vapor generation system, and mechanisms

20 configured to inflate and deflate the balloons as well as control the delivery and removal of integrity gas/fluid and vapor from the device. Additionally, connection lumens 118 can connect device 100 to a vapor generator or vapor source 121, gas/fluid source 122, pressure regulator 124, and flow meter(s) 126. In some embodiments, the vapor generator 121 can include an integrated electronic controller configured to control all aspects of the uterine ablation device,

25 including vapor generation and delivery. In other embodiments, the controller is located in other components of the device, or comprises a separate control unit, such as a computer, smartphone, or tablet that is electronically coupled to the device. Vapor ports 107 near the distal tip 106 of the device can be fluidly coupled to the connection lumens 118 via inflow and outflow lumens (not shown). The vapor ports, inflow and outflow lumens, connection lumens, gas/fluid source,

30 pressure regulator, and flow meters can be configured for testing the integrity of the patient's uterus, proper placement of the device, and verifying the presence of flow between the inflow and outflow lumens of the device. Thermocouples and/or temperature sensors can be disposed in various places within or outside the uterine ablation device, including within or near the inflow and outflow lumens of the device. The thermocouples can be configured to measure temperature

information relating the flow of fluid and/or gas in the device, as well as the temperature within patient body cavities.

[0027] The flow meter can be any flow meter as known in the art, including a thermal mass flow meter, an ultrasonic flow meter, a paddlewheel, or a variable area flow meter. In one embodiment, an ultrasonic flow meter that utilizes transit time and Doppler flow readings is advantageous since it is a non-contact system that does not need to physically interact with the fluid or gas media being employed in the integrity test. An ultrasonic flow meter can be easily adaptable to the exterior dimensions of an inflow lumen. In addition, a drip chamber within the inflow lumen can be used to manually visualize or record drips or flow from the fluid source as the integrity test indicates a sealed uterine cavity. In some uterine procedures, it may be advantageous to use other types of fluid besides saline including Lactated Ringers, non-isotonic solutions for certain electrosurgical procedures, gels, foams, fluids of varying viscosity for some ultrasonographic procedures, or other fluids used in uterine procedures.

[0028] In one embodiment, a one way valve can be placed in the inflow lumen on either side of the flow meter relative to the gas/fluid source. The one way valve can allow for the flow of gas/fluid (e.g., saline) from the gas/fluid source to the device and uterine cavity. The one way valve should not interfere with the operation of the flow meter and its readings. In operation, the uterine cavity is a muscle that can undergo significant contractions during the integrity and patency tests. These contractions can push the fluid retrograde back through the saline lumen and past the flow meter. In doing so, flow meter measurements can become difficult to interpret or may produce sinusoidal waves in the output readings. The placement of the one way valve in the inflow lumen can eliminate retrograde fluid flow and stabilize readings for the flow meter during episodes of uterine contractions.

[0029] Handle 104 can be an ergonomic handle and can include features and controls for using the device (e.g., buttons, levers, indicia for providing feedback for depths of insertion, valves, etc.), including features for controlling inflation of balloons 108, 110, and 112, and for controlling the delivery and removal of integrity test gas/fluid and heated vapor from the device. The handle can also include features and controls for testing the integrity of the patient's uterus, proper placement of the device and verifying the presence of flow between the inflow and outflow lumens of the device.

[0030] The balloons described herein can be any type of flexible balloon, such as rubber, latex, urethane, silicone, PET, LDPE, parylene, nylon, PE, combinations of these polymers, or can be manufactured from any other suitable material as known in the art. It should be noted that in some embodiments, the distal anchor comprises a balloon, but in other embodiments, the distal anchor comprises an expandable anchor or expansion mechanism, such as expandable

frames, filters, nets, or cages, or non-expandable components that increase the diameter of the shaft of the uterine ablation device. For purposes of this disclosure, however, the distal anchor may be referred to as a distal anchor or as a distal balloon.

[0031] Shaft 102 can be configured to deliver a heated vapor from a remote boiler (not shown) through the device and out of vapor ports 107 in distal tip 106. The shaft can also be configured to return vapor that has exited the device, including bodily fluids, uterine materials, and condensate back through the vapor ports and into the shaft. In Fig. 1A, vapor ports 107 are illustrated as including both the vapor delivery and vapor return ports. However, in other embodiments, the vapor delivery ports can be separate and distinct from the vapor return ports. For example, vapor delivery ports are intended to provide an even distribution of heated vapor through a cavity, and may comprise small lumens or holes on the end of the shaft. The vapor return ports, in contrast, are intended to return used vapor and condensate, and may comprise larger slots to prevent blood, tissue, etc. from blocking or clogging the return lumen. The device comprises inflow and outflow gas and/or fluid delivery channels to conduct uterine integrity and patency tests. In some embodiments, the lumens to deliver and return vapor are the same as the channels to deliver and return gas and/or fluid for the uterine integrity and patency tests.

[0032] Referring still to Fig. 1A, uterine ablation device 100 is shown in a collapsed delivery configuration, with distal balloon 108, sealing balloon 110, and positioning balloon 112 deflated to reduce the cross sectional diameter of the device and can be 6 mm in diameter during insertion or smaller. When the device is in the delivery configuration, the reduced profile allows for easier access to through the vagina, cervical canal, and cervix to gain access to the uterus, and provides reduced patient discomfort during insertion. In some embodiments, the outer dimensions of the uterine ablation device are such that introduction of the device into the uterine cavity can be achieved without the need for mechanical or pharmacological dilation of the os prior to device introduction.

[0033] Fig. 1B illustrates the uterine ablation device 100 of Figs. 1A with all three balloons inflated, including distal balloon 108, central sealing balloon 110, and positioning balloon 112. The central balloon can be inflated with a fluid, such as saline, or alternatively, can be inflated with air. Although three balloons are depicted in Fig. 1B, in other variations one, two, four, or more balloons may be provided, and other balloon shapes may be used. The positioning balloon can be inflated with a room temperature medium, a cooled medium, or alternatively, a heated medium. In some embodiments, the central sealing balloon comprises a length along shaft 102 of approximately 15mm to 25mm. The central balloon can be disposed on the shaft between the distal balloon or anchor and the proximal balloon. In some embodiments, the central balloon is adjacent to both the distal balloon and the proximal balloon. In other embodiments, there is a



small gap or space between one or more of the balloons. The length and position of the central balloon on the shaft ensures that when inflated, the central balloon seals the cervix off from the uterus near the internal os, but the balloon does not extend into the uterus or into the vagina of the patient. The central and proximal balloons can comprise any diameter, but preferably should  
5 have a diameter large enough to be able to engage the walls of the cervix and/or the vagina in the average female patient. For instance, the central balloon may have an inflated outer diameter of 10mm and accommodate 9.5 psi of pressure in actual use. The proximal balloon can have a larger diameter, such as 17mm and a lower inflation pressure of 7 psi.

[0034] Fig. 1C is another view of the uterine ablation device 100 including some additional  
10 features. Check valve 144 can connect a vapor delivery conduit to the handle of the device, allowing vapor or other ablation media to be delivered from into the handle. Valve 128c can be a balloon valve configured to control the flow of media through an outflow lumen of the device. The outflow lumen can be used, for example, to remove gas or fluid media from the uterus to a waste container, which will be described in more detail below. Tenaculum stabilizer 146 can  
15 also be seen on the device 100 and can be configured to provide a mechanism to hold a tenaculum and vapor probe handle to maintain position relative to the cervix.

[0035] Fig. 1D is another view of the uterine ablation device 100, further illustrating a tip  
cover 148 disposed over the distal portion of the shaft of the device, including the distal, central, and sealing balloons, and the vapor ports and distal tip. The tip cover can be used in pre-  
20 procedure pressure sensor and balloon checks prior to inserting the device in a patient, which will be described in more detail below.

[0036] One embodiment of a pre-procedure check will now be described. In one  
embodiment, the uterine ablation device is activated by the user and the necessary accessories and device are attached to the vapor generator. Upon attachment of the device handle to the  
25 generator, the accessories are sensed and registered by the electronic controller of the uterine ablation device. Various sensors in the device, such as a water ion sensor and fluid level sensor can provide sensed data to ensure that the proper fluid and amount of fluid have been added to vapor generator and gas/fluid source. The uterine ablation device can then prepare the vapor generator to perform an ablation procedure by heating the generator to its operating levels and  
30 priming steam throughout the internals of the system for a prescribed period of time, ranging from 1 to 5 minutes. The internals can include delivery and return vapor paths within the uterine ablation device. The running of vapor at an appropriate temperature and pressure throughout the internal equipment for a prescribed period of time ensures that the equipment is sterilized for use in a medical procedure. In one embodiment, during this sterilization process, the system can

prevent connection of the vapor delivery conduit of the device to avoid disruption of the sterilization cycle.

[0037] Once sterilization is complete and the vapor delivery conduit is connected, it can be locked into place to prevent inadvertent disruption or disconnection during the rest of the procedure. The vapor delivery conduit of the uterine ablation device can then be pre-heated by allowing vapor to circulate within the vapor delivery conduit. The temperature can then be monitored by a thermocouple in a return path of the generator or a return path of the uterine ablation device. Once completed, the device can notify the user that the generator unit is ready for the next step in the preparation process, or the device can automatically move on to the next step in the preparation process.

[0038] While the sterilization or vapor priming process is occurring, a pressure sensor of the uterine ablation device can be registered, and the pressure reading from the sensor can be compared with a reference sensor within the generator unit. Furthermore, the controller can test the integrity of the three balloons (e.g., distal, central, and sealing) on the device. The system will not allow the process to continue unless all parameters tested meet specifications. In one embodiment, the distal, central, and sealing balloons can be inflated inside the tip cover 148 to a test pressure and checked for leaks. The tip cover can be supplied on the uterine ablation device to also protect the product during shipping and handling within the procedure room. With the tip cover on and the distal, central, and sealing balloons inflated, a chamber created inside the tip cover is pressurized with air through a lumen of the device to compare a pressure sensor reading of the device to a reference pressure sensor in the vapor generator.

[0039] The electronic controller can be configured to prompt the user if the pre-procedural sterilization of the vapor generator or the uterine ablation device is not initiated or completed correctly. The controller can also be configured to prompt the user if the distal, central, and sealing balloons have leaks, or if the pressure and/or temperature sensors indicate a problem with the uterine ablation device or the vapor generator. In some embodiments, the electronic controller can prevent the delivery of condensable vapor to the patient if any of the pre-procedural checks or preparation steps indicate a problem or fail to complete.

[0040] Upon passing the testing of the pre-procedure checks described above, the controller and vapor generator can perform a fluid priming operation to run saline or another fluid through the uterine ablation device in preparation for insertion into the patient. During fluid priming, proper operation of the device can be verified by checking signal strength in the flow sensor and establishing a reference zero flow. The user can then be notified by the system that the uterine ablation device can be inserted into the patient. The tip cover is removed from the distal end of the device prior to insertion.

[0041] Placement of the ablation device of Figs. 1A-1D will now be described. The distal tip of the ablation device can be inserted past an external os into the cervical canal of the patient, and past an internal os of the patient to gain access to the uterus. In one embodiment, the distal balloon can be positioned within the uterus distal to the internal os, the sealing balloon can be positioned at or proximal to the internal os and extending into the cervical canal, and the positioning balloon can be positioned within the cervical canal and extending proximally into or towards the vagina.

[0042] Once the distal tip of the ablation device is disposed within the uterus, just distal to the internal os, the distal balloon can be inflated to the desired pressure. In some embodiments, the balloon can be inflated to a pressure of up to approximately 20 to 30 psi so as to prevent accidental withdrawal of the ablation device from the uterus. It should be noted that at this point, the distal balloon is positioned slightly past the internal os of the cervix. Inflation of the distal balloon can later serve as an anchor to prevent the device from sliding proximally out of the uterus. The user or physician can tug gently on the device to confirm that the distal balloon is in the uterine cavity.

[0043] After inflating the distal balloon, the proximal balloon can be inflated to cause the device to assume a positioned configuration, with the distal balloon fully seated against the internal os and the positioning or proximal balloon expanded within the cervix and extending past the external os into the vagina. As the proximal balloon is inflated, the balloon can expand outwardly from the cervix into the relatively unconstrained space of the vagina, which creates a compression force that pulls the device and the distal balloon proximally to engage against the interior portion of the internal os (also known as the cervical ostium or cervical os).

[0044] Fig. 2 illustrates ablation device 100 inserted into the uterus of a patient with balloons 108, 110, and 112 inflated as described above. Once the distal and proximal balloons are inflated, the central balloon can be inflated to provide additional sealing along the length of the cervical canal. The middle balloon can also be the location of a thermocouple or temperature sensor configured to sense temperature readings in the cervix. If these temperature readings reach a threshold value, such as greater than 44 degrees C, the controller can automatically pause vapor delivery to avoid patient injury. This thermocouple can be tested or verified at the start of the procedure by the controller which queries the thermocouple following placement in the patient to ensure the thermocouple is connected and reading a temperature.

[0045] After positioning the ablation device but prior to delivery of vapor, it can be advantageous to assess the integrity of the uterus to test that the vapor delivery tip of the device is positioned within a sealed uterus and to test that there is flow between the inflow and outflow lumens, by performing an integrity test and a patency test. The amount of fluid and rate in which

it flows into the uterine cavity can provide the physician an indication of the size of the uterine cavity and whether the device is in a false passage. An integrity test can assess that the uterus is sealed, and determine leaks originating from 1) perforations to the uterine wall, or 2) leaks from inadequate sealing at the cervix or 3) leaks from the fallopian tubes.

5 [0046] A second test that made an assessment for patency, referred to as the device lumens patency test or patency test, could provide an indication to the physician whether the device was clogged with debris or placed within a false passage. This additional information to the physician, in conjunction with the integrity test, can provide greater assurance to the physician of device location during “blind” endometrial ablation procedures.

10 [0047] In clinical use, a uterine integrity and patency test could be useful for additional uterine procedures besides uterine ablation procedures such as the implantation of a device, implant, or a diagnostic or therapeutic agent. In these cases, a separate unit or module that can conduct a uterine integrity and patency test, sequentially, separately, or individually, with a separate uterine cavity introducer can be employed without a uterine ablation device or system.

15 [0048] In one embodiment, a uterine integrity test can contain the following elements and steps. Referring to Figs. 1A-1B and Fig. 2, gas/fluid source 122 can be connected to pressure regulator 124 comprising either one regulator or an additional back pressure regulator. The gas/fluid source can contain a gas, such as CO<sub>2</sub>, or inert gases, or a fluid, such as saline, Ringer’s Lactate, non-isotonic solutions, glycerine, and mineral oil for example. The regulator 124 can be  
20 configured to keep the pressure of the external gas source below a safety threshold value. In one embodiment, the safety threshold value can be approximately 70mm Hg. The actual pressure amount or graduation may not be monitored and may not need to be. The fluid or gas from gas/fluid source 122 can be driven at a constant pressure bounded by the safety threshold value (e.g., can be bounded by the maximum pressure the uterus will see during treatment, such as  
25 70mm Hg). In addition, it can be useful to operate a uterine integrity test at a pressure equal to higher than the pressure required for conducting the endometrial ablation or other uterine procedure.

[0049] In one embodiment, gas/fluid pressure can be achieved by elevating the gas/fluid source 122 a height distance above the uterine cavity to create pressure. This height elevation  
30 can be verified by a measuring stick, tape or laser. An example of a clinically used height for a saline bag would be at least 30 inches above the patient’s uterus. At this height, the pressure would be between 50 and 70 mmHg. This pressure is low enough to be below the reported opening pressure of the fallopian tubes. In addition, a pressure sensor within the uterine cavity can verify that the appropriate amount of pressure is being applied for the integrity test and  
35 patency tests. A self-adjusting feedback mechanism can be employed to raise or lower the

pressure of the saline source in response to pressure measurements taken from within the uterine cavity. As an example, this feedback mechanism can raise or lower the height of the saline source in response to the pressure measurements taken from within the uterine cavity.

[0050] In some embodiments, the system can measure a flow rate of gas/fluid exiting the distal lumen of the uterine device or uterine ablation device during the uterine integrity test. This flow rate can also be used to determine the proper pressure or height of the gas/fluid source. For instance, flow rate readings can be taken while the gas/fluid source is at a certain height and the uterine device maintained within a known condition or in free space. As the height of the gas/fluid source is raised or lowered, the flow rate of the gas/fluid will respond accordingly until the gas/fluid source is placed at a height at the desired flow rate, or is pressurized to the desired amount. Likewise, the gas/fluid source can be raised or lowered by a self-adjusting feedback mechanism in response to the measured flow rate.

[0051] In some embodiments, the uterine ablation device can further include a flow meter 126 having a read out mechanism (not shown) to the end user. In one embodiment, the flow meter can comprise an ultrasound sensor, or an optical sensor configured to sense the drip rate of the gas/fluid. In some embodiments, the flow meter can be disposed near distal tip 106 of the device. In other embodiments, the flow meter can be disposed within an outflow lumen of the device. In yet another embodiment, the flow meter can be disposed external to the device but along the flow path between gas/fluid source 122 and the ablation device. The flow meter can be configured to measure and report a flow rate of fluid/gas or vapor as it moves through or exits the uterine ablation device. The read out mechanism can be numerical, graphical, or icon based. Other variations include various audio and visual signals, indicia, qualitative indicia, alarms, and color identifiers. A filter may or may not be attached to the flow meter.

[0052] Referring to Figs. 2 and 3, to perform a uterine integrity test, gas, such as CO<sub>2</sub>, or a fluid, such as saline, can be delivered from the gas/fluid source 122, through a pressure regulator, and through a flow meter 126 into the uterine ablation device 100. As shown in Fig. 3, the gas/fluid can be delivered into the uterus via both inflow lumen 129 and outflow lumen 131. In one specific embodiment, a saline such as 0.9% NaCl can be delivered into the uterus during a uterine integrity test, to determine whether there are leaks in the uterus or cervical canal through which vapor could escape during an ablation procedure. The uterine ablation device 100 can be coupled to an energy generator 124 and controller 123 for uterine ablation therapy. The vapor generator can be, for example, a vapor generator (as shown), but can also be any other type of energy generator, such as an RF energy generator, a cryotherapy generator, etc. Any type of energy modality can be used to ablate or treat the uterus after performing the integrity and patency tests described herein. The energy generator 124 can also include an air pump

configured to supply air to the various balloons in the system, such as the distal, central, and proximal balloons or the balloon valves.

[0053] In one embodiment, a one way valve 127 as seen in Fig. 3 can be located between the flow meter 126 and the uterine ablation device 100. In other variations the one way valve 127 can be located in the handle of the uterine ablation device 100 as well as other components such as the flow meter 126 and valves 128a, 128b, and 128c. The one way valve can reduce or eliminate retrograde flow of saline during uterine contractions. The one way valve is characterized as providing low resistance to flow in one direction (towards the uterine cavity) and high resistance to flow in the retrograde direction (towards the gas/fluid source).

Advantageously the one way valve can stabilize flow values because retrograde flow values are eliminated. By reducing the sinusoidal wave patterns that can be caused by uterine contractions or relaxations, movements by the patient, or inadvertent manipulations of the inflow line or the patient herself by the physician or medical staff, the procedure time is reduced. This filtering out of negative flow values isolates positive components of flow, reduces noise in flow rate values, thereby accelerating the interpretation of flow rate data and reducing procedural time. A fiso airline can connect a reference pressure transducer 125 to the outflow line of the integrity test saline flow pathway.

[0054] A controller of the uterine ablation device, either integrated into the device or into the vapor generator coupled to the device, can be configured to open and close valves 128a, 128b, and 128c to allow gas or fluid to flow from source 122 into the inflow and outflow lumens 129 and 131 of the ablation device 100. Valves 128a, 128b, and 128c can be any type of valve known in the art, such as solenoid valves, inflatable balloons, air cylinders, or electric/hydraulic actuators or cams and gears. During a uterine integrity test, the controller can be configured to open valves 128a and 128b and close valve 128c, to prevent passage of gas/fluid into the waste container 133. This allows gas or fluid to flow from source 122, through flow meter 126, through one way valve 127 and valves 127a and 128b, and into inflow lumen 129 and outflow lumen 131. As the gas or fluid enters the uterus, the flow meter can measure an integrity flow rate of the gas or fluid.

[0055] In one embodiment, the controller of the uterine ablation device or the vapor generator can run an integrity test algorithm to determine if the uterus is sealed. The algorithm can analyze data from the flow meter during the integrity test as gas/fluid is delivered into the uterus. Specifically, the algorithm can analyze a maximum flow rate and a minimum flow rate during an integrity test time window. The integrity test time window can be, for example, a rolling time window of a pre-selected duration. In one specific embodiment, the algorithm analyzes a maximum flow rate and a minimum flow rate continuously during a rolling 15-second

integrity test time window. For each rolling integrity test time window, the minimum and maximum flow rates can be calculated. The difference between the minimum and maximum flow rates in each integrity test time window can be calculated to yield a delta flow value (maximum flow rate minus minimum flow rate), which can be used as an indicator of the stability of flow. For example, the larger the delta flow value, the less stable the flow of gas/fluid, and the smaller the delta flow value, the more stable the flow of gas/fluid. If the maximum flow rate and the delta flow value of gas or fluid stabilizes below an integrity threshold value, the controller can determine that the uterus is sealed. Importantly, the test is comprised of two algorithms that compare flow to an integrity threshold value concurrently with a second algorithm that compares the delta flow value to the integrity threshold value, and uses both of these comparisons to determine the ultimate outcome of the integrity test. The application of both of these comparisons provides greater sensitivity in the test results.

[0056] In some embodiments, this integrity flow rate delta threshold value can be approximately 5 ml/min. Therefore, in some embodiments, a uterus is considered to “pass” the uterine integrity test if both the maximum flow rate and the integrity flow rate delta threshold value are below 5ml/min over a rolling integrity test time window. Alternatively, the test can include different thresholds for maximum flow rate and the delta flow value.

[0057] In some embodiments, the uterine integrity test can run for a pre-set time period. For example, the test can run for 60 seconds, and subsequent rolling 15-second windows can be analyzed to determine if the uterus is sealed during the 60 second time period. In another embodiment, the delta flow value can be defined as a standard deviation of the average flow that is compared to a threshold value. This delta flow value can then be compared to the threshold value to determine if the uterus is sealed.

[0058] In some embodiments, the return channel comprises a valve 128c, such as a solenoid valve, air cylinder, electric/hydraulic actuators, cams and gears or pump/inflatable balloon, which can be activated upon the start of the integrity test to close off the egress of the gas/fluid through the return channel of the uterine ablation device. When the return flow of gas/fluid through the return channel is stopped with the valve, a change of flow can be detected by the flow meter 126 on the input line. In addition to determining if there is a leak or if the device is positioned properly, the specifics of the changes in flow (e.g., how the flow reacts to closing of the return line with the valve) can provide the following the indications in some cases: a) the size of the uterine cavity; and b) the presence of a leak or lack of integrity in the system. For instance in clinical use with uteri of varying sizes, an integration under the graphical curve of flow rate versus time provides a volume assessment of the size of uterine cavity. The amount of volume can provide the physician information not only on the size of the uterus, but whether the

device is improperly embedded in a false passage (smaller volume amount) or in the peritoneal cavity (larger volume amount).

[0059] Immediately after performing the integrity test above, the amount of flow in the inflow and outflow channels can be measured in a patency test and used to determine the presence of an obstruction that may affect the flow of vapor during the ablation procedure.

Based on this determination or patency test, the device may be repositioned or replaced prior to delivery of vapor. For example, in one embodiment, referring to Fig. 4, a method of performing a patency test can comprise delivering gas or fluid from inflow lumen 129 of the uterine device into the uterus, also referred to as the fluid infusion tip, removing gas or fluid from the uterus with outflow lumen 131 of the uterine device, also referred to as the fluid outflow tip, and determining that the uterine device is not clogged or embedded in tissue if a flow rate of gas or fluid is observed in the flow meter of the inflow lumen of the uterine device. In Figs. 3-4, valves 128a and 128b control the flow of gas/fluid to the uterine ablation device 100 and valve 128c control the flow of gas/fluid from the outflow lumen 131 into the outflow canister or waste container 133. Control of the valves 128a and 128b and 128c can be performed by a separate controller and software unit shown as 123.

[0060] If it has been determined that the uterus is sealed based on the integrity test performed and described in Fig. 3, the controller can also be configured to perform a patency test. In one embodiment, referring to Fig. 4, the controller can be configured to open valves 128b and 128c, but close valve 128a. This allows gas or fluid to flow from source 122, through flow meter 126, through one way valve 127 and valve 128b, and into inflow lumen 129. Gas or fluid can be removed through outflow lumen 131, through valve 128c, and into a waste container 133. As the gas or fluid enters and is removed from the uterus, the flow meter can measure a patency flow rate of the gas or fluid. If the patency flow rate is maintained above a patency flow rate threshold value, the controller can determine that the device is not clogged or embedded into tissue. In some embodiments, observing or measuring a flow of fluid or gas in outflow lumen 131 can be used to determine that the device is not clogged or embedded in tissue. A flow rate above a patency test threshold during a rolling patency test time window can indicate that the lumens are not clogged or that the distal end of the uterine ablation device is not embedded into tissue.

[0061] In one specific embodiment, the patency test threshold can be greater than 5 ml/min, and the rolling patency test time window can be a 5 second time period. Thus, the flow meter can measure the patency flow rate in rolling patency time windows (e.g., rolling 5 second periods) and the controller can analyze the measured rate. If the patency flow rate is maintained above the patency test threshold (e.g., 5 ml/min) during a rolling patency time window, then the



patency test is considered passed and the test can be stopped. Passing the patency test indicates that the uterine ablation device is not obstructed or placed in false passage. If the patency test threshold is not satisfied, the physician should repeat the insertion steps and/or repeat the integrity test and patency test prior to initiating uterine ablation. When the patency flow rate is below the threshold of 5 ml/min during the rolling patency test time window, the uterine ablation device may need to be repositioned.

[0062] During the transition from the end of integrity test to the start of the patency test, the uterine cavity can be substantially filled with the gas/fluid provided during the integrity test. As described above, the closed outflow valve during the integrity test prevents gas or fluid from exiting the uterine cavity into the waste container 133. In one embodiment, it is desirable for valve 128c to be opened only partially in the range of 20-50% open for a flow rate greater than 5 ml/min and less than 40 ml/min so the uterine cavity distension achieved during the integrity test is temporarily maintained when the patency test checks for open flow through the uterine ablation device. Certain types of valves are better suited for partial opening. For example, balloon valves can be pulsed at various duty cycles to partially open the valve. The higher the duty cycle, the more quickly the valve can be opened. Partial opening of the valve prevents the uterine cavity from collapsing too quickly around the tip of the uterine ablation device which, in some instances, may cause a false positive failure of the patency test. In one embodiment, partial opening of the valve can be achieved by pulsing the opening of the valve at a specified duty cycle until flow through the vapor probe begins, or alternatively until the uterine pressure begins to drop. In another embodiment, the valve 128c can be opened rapidly just until flow through the valve begins. This rapid drop opening of the valve can be achieved by pulsing the valve initially with a high duty cycle, then shortening the pulsing (or lowering the duty cycle) as the valve approaches the range where flow through the valve begins. Once the patency flow rate increases above a threshold (or by a specific rate of increase), the valve can be maintained.

[0063] In one specific embodiment, the valve 128c can be a balloon filled to as much as 20 psig to occlude the tubing leading to waste container 133. The balloon valve can be pulsed open for up to 40 msec every 200 msec until the balloon pressure falls to as low as 5 psig. The valve opening time can then be reduced even further until the balloon pressure falls to between 3-4 psig. The valve can continue to be pulsed until flow increases to a level of 0.20 ml/min or until flow rises above the threshold value (e.g., above 5 ml/min).

[0064] Fig. 5 shows one specific embodiment where the outflow valve (e.g., valve 128c from Fig. 3) comprises an inflatable balloon. Inflating the balloon causes obstruction of the outflow lumen of the uterine ablation device, and deflating the balloon allows flow out of the device and into the waste container. Fig. 5 shows the outflow valve pressure 136 (mmHg), the flow 138

(ml/min) through the uterine ablation device, and the uterine cavity pressure 140 (mmHg) during a typical patency test. At the completion of a typical integrity test, the cavity pressure should be above 52mmHg, the outflow valve pressure about 20psi, and the flow through the uterine ablation device close to zero. In the first few seconds of the patency test, the outflow valve pressure decreases rapidly, then the rate of deflation decreases to zero as flow through the device begins, shown by arrow 142. The cavity pressure drops gradually as flow increases. The patency algorithm can run concurrently with deflation of the outflow valve. In some embodiments, the deflation period of the outflow valve is typically from 3 to 40 msec.

[0065] Upon successful completion of the integrity and patency test, the controller of the uterine ablation device can automatically begin the treatment process. The treatment process can comprise delivering vapor through the outflow lumen and vapor ports of the uterine ablation device. To collect condensate, debris, and byproducts of the ablation process, and free flowing vapor, the return lumen can include a micro-porous filter at its distal end. A fiber optic pressure sensor can be located within the filter to ensure that the treatment pressure stays within a prescribed range. Pressure control can be achieved through a series of valves that automatically engage and disengage to maintain the treatment pressure within the prescribed range. If the pressure sensor records values too low or too high, the device can automatically pause vapor delivery and alert the physician of the condition that caused discontinuation of vapor delivery.

[0066] The uterine ablation device can further include an interrupt button and on-screen button to pause the procedure, if necessary, at any time during the procedure. As described above, a thermocouple can be located in the lumen that exhausts vapor from the uterine cavity. This thermocouple can provide post-procedural feedback on the temperature of the byproducts flowing through the outflow line and into the waste container. This temperature provides an indication of the amount of vapor that has circulated through the cavity and exited through the return lumen and outflow tubing.

[0067] Ablation is achieved by imposing elevated temperatures on the endometrial tissue for long enough duration to bring about cell necrosis. Water vapor's unique feature is the high concentration of thermal energy that is characteristic of the vapor phase. The easy extraction of the concentrated thermal energy by means of condensation enables high rates of heat transfer to be applied to the walls of the uterine cavity. The vapor is delivered in the uterine cavity at a prescribed pressure range and at a saturation temperature that is approximately one degree above the normal boiling temperature of water. The uterine surface temperature is below that of the vapor, giving rise to sustained condensation.

[0068] The rate of utilized vapor flow can be auto-regulated by the rate of condensation and the rate at which heat penetrates the uterine tissue. If the rate of delivered vapor flow exceeds

the auto-regulated value, it will leave the uterine cavity through the return lumen and outflow tubing of the device. The thermal energy deposited on the uterine wall by condensation penetrates the uterine tissue by two transport mechanisms. One of these is heat conduction, which is a molecular-level process that is indigenous to all media. A second transport

5 mechanism is due to blood perfusion, which is a seepage-like process. Heat is carried by the perfusing blood, also called the "heat sink" effect. Therapeutic cell necrosis is achieved both maintaining elevated temperatures at the inner surface of the uterus and a long enough period to achieve the desired depth of ablation. Since an excess of vapor at a prescribed pressure and period of time is supplied during the procedure to achieve ablation, a large range of cavity sizes  
10 can be treated uniformly to the desired depth without the need to adjust the treatment parameters. Uniformity of coverage within the cavity is achieved by the continuous supply of steam that circulates within the cavity prior to exiting through the return lumen and outflow tubing.

[0069] At the completion of treatment at a predetermined time, which can be less than 140 seconds in one embodiment, the delivery of vapor is stopped by the controller and generator unit.  
15 After the cessation of vapor delivery, the generator can query the intrauterine pressure sensor to ensure that internal pressure values have now dropped. Upon reaching a satisfactory level, the system automatically deflates the balloons holding the device in the cervix and the physician can be notified that the delivery device is now ready for removal from the patient. Following this step, the system can unlock the uterine ablation device from the generator and the delivery  
20 device can be discarded.

[0070] Once treatment is completed, the system also prompts the end user if another patient is about to be treated. If so, the system repeats the steps described above. If there is not another patient to be treated, the excess fluid can be drained from the generator.

[0071] As for additional details pertinent to the present invention, materials and  
25 manufacturing techniques may be employed as within the level of those with skill in the relevant art. The same may hold true with respect to method-based aspects of the invention in terms of additional acts commonly or logically employed. Also, it is contemplated that any optional feature of the inventive variations described may be set forth and claimed independently, or in combination with any one or more of the features described herein. Likewise, reference to a  
30 singular item includes the possibility that there are plural of the same items present. More specifically, as used herein and in the appended claims, the singular forms "a," "and," "said," and "the" include plural referents unless the context clearly dictates otherwise. It is further noted that the claims may be drafted to exclude any optional element. As such, this statement is intended to serve as antecedent basis for use of such exclusive terminology as "solely," "only" and the like in  
35 connection with the recitation of claim elements, or use of a "negative" limitation. Unless

defined otherwise herein, all technical and scientific terms used herein have the same meaning as commonly understood by one of ordinary skill in the art to which this invention belongs. The breadth of the present invention is not to be limited by the subject specification, but rather only by the plain meaning of the claim terms employed.

## CLAIMS

## WHAT IS CLAIMED IS:

- 5 1. An endometrial vapor ablation system, comprising:  
a vapor generator;  
a uterine ablation device coupled to the vapor generator, the uterine ablation device  
configured for insertion into a uterine cavity of a patient;  
a control system coupled to the uterine ablation device and the vapor generator, the  
10 control system being configured to automatically perform a pre-procedural sterilization of the  
vapor generator and the uterine ablation device, the control system configured to prep the uterine  
ablation device for vapor delivery, the control system configured to test safety components of the  
uterine ablation device prior to and after insertion into the patient, the control system configured  
to monitor the application of vapor energy, monitor temperature reading at identified locations  
15 on the uterine ablation device, and the control system being configured to implement shutdown  
of the vapor generator following vapor delivery.
2. The system of claim 1, wherein the electronic controller is configured to prompt an end  
user of the uterine ablation device if the pre-procedural sterilization fails to initiate or complete.  
20
3. The system of claim 1, wherein the electronic controller is configured to prompt an end  
user of the uterine ablation device if the testing of the safety components indicates an error in the  
system.
- 25 4. The system of claim 1, further comprising a distal anchor balloon, a central sealing  
balloon, and a proximal positioning balloon disposed on a shaft of the uterine ablation device.
5. The system of claim 4, further comprising a temperature sensor positioned near the  
central sealing balloon and configured to measure a temperature inside an endocervix of the  
30 patient.
6. The system of claim 4, further comprising one or more pressure sensors configured to  
measure a pressure of each of the distal anchor balloon, the central sealing balloon, and the  
proximal positioning balloon.  
35

7. A method of providing vapor therapy to a uterus of a patient, comprising the steps of:  
sensing a fluid level of a vapor generator;  
heating the vapor generator to prepare a condensable vapor;  
priming the condensable vapor through a uterine ablation device coupled to the vapor  
5 generator to sterilize the uterine ablation device;  
comparing a pressure sensor value within the uterine ablation device to a pressure sensor  
value within the vapor generator;  
verifying the integrity of a plurality of positioning and sealing balloons of the uterine  
ablation device;  
10 delivering condensable vapor from the uterine ablation device into the uterus; and  
during the delivering step, monitoring a temperature in a cervix of the patient and  
automatically terminating delivery of condensable vapor if the temperature exceeds a threshold  
value.
- 15 8. The method of claim 7, wherein the threshold value comprises 44 degrees C.





Fig. 1A








Fig. 1C

Fig. 1D





Fig. 2






Fig. 4

**PATENCY TEST  
SALINE FLOW PATHWAY**





Fig. 5



## INTERNATIONAL SEARCH REPORT

International application No.

PCT/US15/32009

<b>A. CLASSIFICATION OF SUBJECT MATTER</b> IPC(8) - A61B 17/42, 18/04 CPC - A61B 17/42, 18/04, 2018/00005, 2018/00011, 2018/00017 According to International Patent Classification (IPC) or to both national classification and IPC		
<b>B. FIELDS SEARCHED</b> Minimum documentation searched (classification system followed by classification symbols) IPC(8): A61B 17/42, 18/00, 18/04 CPC: A61B 17/42, 18/04, 2018/00005, 2018/00011, 2018/00017 Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched Electronic data base consulted during the international search (name of data base and, where practicable, search terms used) PatSeer (US, EP, WO, JP, DE, GB, CN, FR, KR, ES, AU, IN, CA, INPADOC Data); Google/Google Scholar; PubMed/MEDLINE; IP.com: vapor*, ablat*, uterus*, uterine*, peak*, threshold*, endometr*, vagina*, cervix*, temperature*.		
<b>C. DOCUMENTS CONSIDERED TO BE RELEVANT</b>		
Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	US 2009/0149846 A1 (HOEY, M et al.) June 11, 2009; abstract; figures 3, 7-8, 11, 18a; paragraphs [0009], [0019], [0060], [0081], [0084], [0086], [0092], [0097], [0109], [0111], [0113], [0119], [0124]-[0126], [0127], [0134], [0138]	1-3 -- 4-8
Y	US 2013/0090572 A1 (PELIKS, RB et al.) April 11, 2013; abstract; figures 1B-2; paragraphs [0028], [0051], [0058], [0060], [0072], [0074], [0092]; claim 1	4-8
Y	US 2005/0171574 A1 (RUBINSKY, B et al.) August 04, 2005; abstract; paragraphs [0020], [0102]; claims 17-18	8
A	US 5277201 A (STERN, RA) January 11, 1994; entire document	8
<input type="checkbox"/> Further documents are listed in the continuation of Box C. <input type="checkbox"/> See patent family annex.		
* Special categories of cited documents: "A" document defining the general state of the art which is not considered to be of particular relevance "E" earlier application or patent but published on or after the international filing date "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified) "O" document referring to an oral disclosure, use, exhibition or other means "P" document published prior to the international filing date but later than the priority date claimed "T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention "X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone "Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art "&" document member of the same patent family		
Date of the actual completion of the international search 07 AUGUST 2015 (07.08.2015)		Date of mailing of the international search report <b>25 AUG 2015</b>
Name and mailing address of the ISA/ Mail Stop PCT, Attn: ISA/US, Commissioner for Patents P.O. Box 1450, Alexandria, Virginia 22313-1450 Facsimile No. 571-273-8300		Authorized officer Shane Thomas PCT Helpdesk: 571-272-4300 PCT OSP: 571-272-7774