#### (12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

# (19) World Intellectual Property Organization

International Bureau

# (43) International Publication Date 1 March 2012 (01.03.2012)



## (10) International Publication Number WO 2012/027490 A2

(51) International Patent Classification: A61F 2/44 (2006.01)

(21) International Application Number:

PCT/US2011/048992

**English** 

(22) International Filing Date:

24 August 2011 (24.08.2011)

(25) Filing Language:

(26) Publication Language: English

(30) Priority Data:

61/376,626 24 August 2010 (24.08.2010) US 61/526,630 23 August 2011 (23.08.2011) US

- (71) Applicant (for all designated States except US): STOUT MEDICAL GROUP, L.P. [US/US]; 410 East Walnut Street, Suite 8, Perkasie, Pennsylvania 18944 (US).
- (72) Inventors; and
- (75) Inventors/Applicants (for US only): GREENHALGH, E. Skott [US/US]; 820 Penllyn Pike, Lower Gwynedd, Pennsylvania 19002 (US). ROMANO, John-Paul [US/US]; 59 Skyline Drive, Chalfont, Pennsylvania 18914 (US).
- Agents: LEVINE, David A. et al.; Levine Bagade Han LLP, 2400 Geng Rd, Suite 120, Palo Alto, California 94303 (US).
- (81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM,

AO, AT, AU, AZ, BA, BB, BG, BH, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IS, JP, KE, KG, KM, KN, KP, KR, KZ, LA, LC, LK, LR, LS, LT, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PE, PG, PH, PL, PT, QA, RO, RS, RU, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.

(84) Designated States (unless otherwise indicated, for every kind of regional protection available): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).


#### **Declarations under Rule 4.17:**

as to applicant's entitlement to apply for and be granted a patent (Rule 4.17(ii))

#### Published:

without international search report and to be republished upon receipt of that report (Rule 48.2(g))

(54) Title: SUPPORT DEVICE AND METHOD FOR USE



(57) Abstract: Devices and methods for orthopedic support are disclosed. The device can have a first rigid section hingedly attached to a second rigid section. The device can be curved or rotated around obstructions along an access path to a target site. The device can be delivered to an intervertebral location in a patient.





TITLE OF THE INVENTION 1 SUPPORT DEVICE AND METHOD FOR USE 2 E. Skott Greenhalgh 3 4 John-Paul Romano 5 CROSS-REFERENCE TO RELATED APPLICATIONS 6 [0001] This application claims priority to U.S. Provisional Patent Application No. 61/376,626, 7 filed 24 August 2010, and to U.S. Provisional Patent Application No. 61/526,630, filed 23 8 August 2011, both of which are incorporated by reference herein in their entireties. 9 10 11 BACKGROUND OF THE INVENTION 1. Field of the Invention 12 [0002] A device, such as a flexible spinal fusion cage, which can articulate (bend) in such a way 13 that it will be able to be implanted from a lateral approach into L4-L5 and L5-S1 is disclosed. 14 15 2. Description of the Related Art 16 [0003] Typical lateral approach fusion implants (e.g., Discover XLIF, by NuVasive, Inc., San 17 Diego, CA; and the Direct Lateral Interbody Fusion (DLIF) by Medtronic, Inc., Minneapolis, 18 MN) are not able to implant their fusion cages for two reasons. 19 [0004] First, boney obstacles can impair access. Figures 1a and 1b illustrate the pelvis and lower 20 spine including the Ilium 2, sacrum S1, and lower lumbar vertebrae L3, L4 and L5. Figures 1a 21 and 1b show the challenge of gaining lateral access to the L4-L5 and the L5-S1 intervertebral 22 spaces. The position of the Hium 2 obstructs the direct lateral access pathway. 23 [0005] Figure 2 illustrates windows 4a and 4b or channels which some doctors create during 24 implantation. The windows 4a and 4b are created through the Ilium to gain direct line of site 25 access to the L4-L5 and L5-S1 intervertebral spaces, respectively. This is a highly invasive 26 approach, creates significant tissue damage, particularly to the Hium and surrounding soft tissue, 27 and requires significant surgical skill. 28 [0006] Second, the steep approach angle (8a for the L4-L5 intervertebral space and 8b for the 29 30 L5-S1 intervertebral space), as measured from a transverse plane along the approach path (10a for the L4-L5 intervertebral space and 10b for the L5-S1 intervertebral space) of a tissue 31

1

retractor relative to the location of the fusion site, can cause problems, as illustrated in Figures 3 1 and 4. The approach paths 10a and 10b pass through the skin surface 12. The tissue retractor 2 3 used in lateral fusion surgery provides line of site access to the disk space requiring a fusion cage insertion. The tissue retractor holds tissue out of the way of the procedure. The tissue retractor is 4 also used to create a working channel to pass tools through, protect neural tissue, and anchor to 5 the superior and inferior vertebral bodies relative the disk space requiring fusion. The volume 6 within the pelvis and inferior to the dashed demarcation line 6 along a transverse plane is very 7 hard if not impossible to reach with a direct lateral approach due to the Ilium. Even if the 8 retractors are tilted as shown by the demarcation line 6, the ability to insert an implant that is the 9 length of the end plates of the L4 or L5 vertebral bodies would be very difficult due to 10 obstruction of the Ilium among other factors. 11 [0007] Furthermore, with the retractor positioned along the approach path 10a or 10b plane and 12 angled direction, the angle formed between the retractor and the vertebral body end plates would 13 make inserting a monolithic, inflexible fusion cage 14 or implant into the L5-S1 intervertebral 14 space difficult if not virtually impossible due to obstruction of the surrounding hard and soft 15 16 tissue, as illustrated by Figure 5a. A typical lateral fusion cage or implant width 16 is the width of the end plate 18 along the adjacent disk. The implant 14 can not turn the corner at the pivot 17 point 20 at the lateral and/or anterior edge of the L5-S1 intervertebral space. 18

19

20

#### SUMMARY OF THE INVENTION

[0008] Support or fixation devices and methods for access, controlling (e.g., steering or rotating, 21 and driving or translating) implants, and modifying the configuration of implants are disclosed. 22 The device can be an implantable fixation device, such as a flexible and/or articulatable fusion 23 cage. The device can articulate and/or bend so the device can make the turn around the L5-S1 24 intervertebral space. The implant can flex and/or articulate. For example, the implant can have 25 hinges and/or be flexible (e.g., have significantly elastic structural components). 26 [0009] Articulation tools are disclosed that can be used to implant the device. The articulation 27 tools can articulate the device and/or allow the device to articulate. For example, the connection 28 between the articulation tool and the implant can bend, flex, steer, or combinations thereof. The 29 articulation tools can be used to debride or clear out the disk space. 30

1 [0010] An oblique curved access tool or device can be used. The device can be delivered to the

- 2 intervertebral space along an oblique approach path, not perpendicular to the spine. The oblique
- approach can provide an access path from lateral skin to the L5-S1 disk space, and can curve
- 4 tangent to the Hium. A large working channel through the soft tissue can be created. The
- 5 oblique access tool can move soft tissue out of the way to create the working channel. The
- 6 oblique approach can reduce the access-tool-to-disk-space approach angle.
- 7 [0011] A biological implant support device for providing orthopedic support is disclosed. The
- 8 device can be articulatable or flexible. The device can have a first rigid section at a first terminal
- 9 end of the device. The first rigid section can have a first top plate and a first bottom plate. The
- device can have a second rigid section having a second top plate and a second bottom plate. The
- 11 first rigid section can be rotatably attached to the second rigid section. The top and bottom plates
- can be configured to interface with hard tissue.
- 13 [0012] A method for inserting a support device to a target site in a spine adjacent to a first
- vertebra is disclosed. The method can include creating a channel through a non-vertebral bone.
- 15 The method can include inserting a first rigid section of the device through the channel and into
- the target site. The method can include inserting a second rigid section of the device through the
- channel. The method can include rotating the second rigid section of the implant with respect to
- the first rigid section. The first rigid section can be hingedly attached to the second rigid section.
- 19 The method can include inserting the second rigid section of the implant into the target site.
- 20 [0013] Creating the channel can include drilling the tissue with a flexible drill. The non-
- 21 vertebral bone can be the pelvis, such as the ilium and/or the sacrum.
- 22 [0014] A method for inserting an implant to a target site between a first vertebra and a second
- 23 vertebra is disclosed. The method can include creating a first channel through the ilium. The
- 24 method can include creating a second channel through the sacrum. The first channel can be
- 25 aligned with the second channel. The method can include inserting a first rigid section of the
- 26 implant through the first channel and the second channel into the target site. The method can
- 27 include rotating a second rigid section of the implant with respect to the first rigid section,
- 28 wherein the first rigid section is hingedly attached to the second rigid section. The method can
- 29 include inserting the second rigid section of the implant into the target site. The second channel
- 30 can pass through a port formed in vertebral endplate. The device can be inserted through the port
- 31 in the vertebral endplate and articulate as the device is delivered into the target site.

1 2

#### BRIEF DESCRIPTION OF THE DRAWINGS

- 3 [0015] Figures 1a and 1b are anterior and lateral views, respectively, of the lower lumbar and
- 4 sacral spine and pelvis with the Ilium shown in phantom lines in Figure 1b.
- 5 [0016] Figure 2 is a lateral view of the lower lumbar spine with windows cut through the Ilium.
- 6 [0017] Figures 3 and 4 are anterior and lateral views, respectively, of the lower spine and pelvis
- 7 along with approach paths into the intervertebral spaces.
- 8 [0018] Figure 5a is an anterior close-up view of the lower spine and pelvis with an approach of a
- 9 monolithic implant.
- 10 [0019] Figure 5b illustrates a variation of the implantable device.
- 11 [0020] Figures 5c and 5d illustrate a variation of a method of delivering the device of Figure 5b
- into the L5-S1 space.
- 13 [0021] Figures 6 through 8 are anterior, perspective and lateral views, respectively, of a variation
- of the approach path for delivering the implant into the intervertebral space.
- 15 [0022] Figures 9a through 9d illustrate variations of the device in various configurations. An x-
- axis, y-axis and z-axis are also shown for orientation with the x-axis disposed along the
- 17 longitudinal axis of the device.
- 18 [0023] Figures 10a and 10b illustrate various configurations of a variation of the device in a
- 19 steering tube with the tube shown as see-through for illustrative purposes.
- 20 [0024] Figures 10c through 10e illustrate various configurations of a variation of the device on
- 21 steering rails attached to the lateral outside of the device.
- 22 [0025] Figures 11a through 11c illustrate various configurations of a variation of the device on a
- 23 steering rail attached to the inside of the device.
- 24 [0026] Figures 12a through 12f are cross-sections of various steering rails, or along the length of
- 25 the same steering rail.
- 26 [0027] Figure 13 illustrates a method for deploying the device into the L5-S1 intervertebral
- 27 space.
- 28 [0028] Figures 14a and 14b illustrate various configurations of a variation of the device in a
- 29 steering slide.
- 30 [0029] Figures 15a and 15b are top and side views of a variation of the device with parallel
- 31 hinges.

4

- 1 [0030] Figure 16 is a top view of a variation of the device with non-parallel hinges.
- 2 [0031] Figures 17a through 17f are side views of variations of the device.
- 3 [0032] Figures 18 and 19 are perspective views showing the orientation of variations of living
- 4 hinges within devices.
- 5 [0033] Figures 20a through 20c are perspective, top and front views, respectively, of a variation
- 6 of the device in a straight or flat configuration.
- 7 [0034] Figures 21a through 21c are perspective, top and front views, respectively, of the device
- 8 of Figures 20a through 20c in an articulated configuration.
- 9 [0035] Figures 22a through 22c are perspective, top and front views, respectively, of a variation
- of the device in a straight or flat configuration.
- 11 [0036] Figures 23a through 23c are perspective, top and front views, respectively, of the device
- of Figures 22a through 22c in an articulated configuration.
- 13 [0037] Figure 24 illustrates the lower spine and pelvis.
- 14 [0038] Figures 25 through 28 illustrate a variation of a method of delivering the device to a
- 15 target site.
- 16 [0039] Figures 29 through 31 illustrate views through the transverse plane from a superior
- 17 location, the sagittal plane from a lateral location, and the coronal plane from an anterior
- 18 location, respectively, of a variation of the location of the transosseous delivery channel.
- 19 [0040] Figures 32a through 32d illustrate a superior view of a variation of a method of delivering
- 20 the device showing the iliac and sacrum, but not the L5-S1 disc or remainder of the spine for
- 21 illustrative purposes.
- 22 [0041] Figures 33a through 33d illustrate a posterior perspective view of a variation of a method
- 23 of delivering the device showing the iliac and sacrum, but not the L5-S1 disc or remainder of the
- 24 spine for illustrative purposes.

25 26

#### DETAILED DESCRIPTION

- 27 [0042] Support or fixation devices and methods for access, controlling (steering) implants, and
- 28 modifying implants are disclosed. The support device disclosed herein can be used to treat one
- or more osseous structures in the body including the L4-L5 and L5-S1 region of the spine. The
- 30 device can be used with known methods of accessing the vertebrae of the spine such as the L4-
- 31 L5 and L5-S1 regions with posterior, anterior, or lateral approaches, or combinations thereof.

1 [0043] The device can be an implantable fixation device, such as a flexible fusion cage. The

- 2 device can be delivered into an intervertebral space, for example, to provide structural support
- 3 between the adjacent vertebrae. The device can fuse the vertebra adjacent to the specific
- 4 intervertebral space. A discectomy can be performed at the target implant site before or during
- 5 delivery of the implant.
- 6 [0044] Figures 5a through 5c illustrate that the device can be articulatable or flexible. The
- 7 implantable device 10 can be used to support and/or fix structures between adjacent vertebrae,
- 8 such as between the L4 and L5 vertebrae or between the L5 and S1 vertebrae. The implantable
- 9 device 10 can be articulatable and/or flexible so as to navigate sharp anatomical turns, such as
- the L4-L5 or L5-S1 intervertebral space. The implantable device 10 can be rigidly lockable or
- 11 can remain flexible or articulatable at all times. The implantable device 10 can be rigidly locked
- for example using a delivery tool, e.g., wires, sheaths, guides, or combinations thereof, for
- example, for additional stability. Such surgical delivery tools, alone or in combination, may add
- 14 axial strength and stability before during or after pressing the implantable device 10 into the
- 15 targeted intervertebral disc space.
- 16 [0045] Figure 5b illustrates that the implantable device 14 can have first, second, third, and
- fourth segments 22a through 22d. Each of the segments 22a, 22b, 22c, and 22d can be attached
- to the adjacent segment at a flex point or articulatable hinge 24a, 24b, and 24c, respectively. The
- device 14 can articulate and/or bend at the hinges 24.
- 20 [0046] Figures Sc and Sd illustrate that the device 14 can be delivered into the L5-S1
- 21 intervertebral space. The device 14 can make the turn around the L5-S1 intervertebral space,
- such as at the pivot point 20, by articulating or flexing.
- 23 [0047] Figures 6 through 8 shows illustrate a curved implant pathway or approach path 10c. An
- 24 articulation tool can be used to push (e.g., impact), pull, control or combinations thereof, the
- 25 implant 14. The implant 14 can articulate and/or flex during delivery. The implant can have
- single or multiple hinges, a flexible shaft, laser slots (e.g., in a tube to act as hinges) or
- 27 combinations thereof.
- 28 [0048] The approach path 10c can be tangential to the medial surface of the Ilium along a
- 29 portion of the length of the approach path 10c. A portion of the length of the approach path 10c
- 30 can be linear and a portion of the length of the approach path 10c can be curved. The entire
- 31 approach path 10c can be linear or curved. A portion of the length of the approach path 10c can

track (i.e., follow the same shape of) the medial surface of the Ilium. The approach path 10c can

- 2 contact the medial surface of the Ilium 2. The approach path 10c can be non-perpendicular or
- 3 perpendicular to the longitudinal axis 27 of the spine where the approach path 10c enters the
- 4 intervertebral space L4-L5 or L5-S1.
- 5 [0049] The approach-Ilium gap 26 can be measured between the approach path 10c and the
- 6 closest medial surface of the Ilium 2. The approach-Ilium gap 26 can be perpendicular to the
- 7 approach path 10c and the Ilium 2, for example when the approach path 10c is tracking the
- 8 medial surface of the Ilium 2. The approach-Ilium gap 26 can be from about 0 mm to about 15
- 9 mm along the length of the approach path 10c where the approach path is tracking the medial
- surface of the Ilium 2, more narrowly from about 0 mm to about 10 mm, yet more narrowly from
- about 2 mm to about 8 mm.
- 12 [0050] The approach path 10c can be curved in all three dimensions (e.g., in the transverse
- plane, sagittal plane and coronal plane), or any combination thereof and straight in the remaining
- 14 dimensions.
- 15 [0051] Figure 9a through 9d illustrate that variations of hinges 24a and 24b between the
- segments 22a, 22b and 22c can allow the implant 14 to articulate. The implant 14 can have
- controlled angulation or articulation (i.e., with discrete, defined built-in stopping points or stops)
- or free angulation or articulation (i.e., with no stops).
- 19 [0052] Figure 9a illustrates that the hinges 24a and 24b can be oriented in parallel with the z-
- 20 axis. The hinges can have a single degree of rotational freedom. The segments 24, 24b and 24c
- 21 can articulate by rotating about the z-axis with respect to each other. The hinges 24a and 24b
- can be near the top (as shown), near the bottom, in the middle with respect to the y-axis, or
- 23 combinations thereof of the device 14.
- 24 [0053] Figure 9b illustrates that the hinges 24a and 24b can be oriented in parallel with the x-
- 25 axis. The segments 24, 24b and 24c can articulate by rotating about the x-axis with respect to
- each other. The hinges 24a and 24b can be near the front (as shown), near the rear, in the middle
- with respect to the z-axis, or combinations thereof of the device 14.
- 28 [0054] Figure 9c illustrates that the hinges 24a and 24b can be oriented in parallel with the y-
- 29 axis. The segments 24, 24b and 24c can articulate by rotating about the y-axis with respect to
- 30 each other. The hinges 24a and 24b can be near the front (as shown), near the rear, in the middle
- 31 with respect to the z-axis, or combinations thereof of the device 14.

- 1 [0055] Figure 9d illustrates that the hinges 24a and 24b can be ball-in-socket hinges allowing
- 2 three rotational degrees of freedom, or a combination of the three hinges described in Figures 9a
- 3 through 9c, allowing two or three degrees of freedom. The segments 24, 24b and 24c can
- 4 articulate by rotating about the x-axis, and/or y-axis, and/or z-axis with respect to each other.
- 5 The hinges 24a and 24b can be near the front (as shown), near the rear, in the middle with
- 6 respect to the z-axis, near the top, near the bottom, in the middle with respect to the y-axis (as
- 7 shown), or combinations thereof of the device 14.
- 8 [0056] The first hinge 24a can be located in a different location and/or with a different than the
- 9 second hinge 24b. For example, the first hinge 24a can be oriented in parallel with the z-axis,
- allow rotation about the z-axis and be located near the top of the device 14, and the second hinge
- 24b can be oriented in parallel with the x-axis, allow rotation about the x-axis, and be located
- near the middle of the device 14 with respect to the z-axis.
- 13 [0057] Figures 10a and 10b illustrate that the device 14 can have an outer steering sheath or tube
- 14 28. The device 14 can be fixed to the steering tube 28 or can slide along the steering tube 28.
- 15 The steering tube 28 can be articulatable and/or flexible, as shown by the arrow in Figure 10b
- and the various configurations of the tube 28 between Figures 10a and 10b. The articulation or
- 17 flexion of the steering tube 28 can be controlled, for example by delivering controlled tension to
- tensile control wires in the walls of the steering tube 28.
- 19 [0058] The steering tube 28 can be positioned at the target deployment site. For example, the
- steering tube 28 can be placed in the intervertebral space and can remain in the intervertebral
- 21 space post-surgery, or the steering tube 28 can be removed from the intervertebral space and the
- device 14 can be deployed from the tube 28 and the device 14 can be left in the intervertebral
- 23 space.
- 24 [0059] Also for example, the distal end of the steering tube 28 can be positioned at the entrance
- 25 to the intervertebral space and/or rested on the inferior and/or superior vertebral body end plate
- adjacent to the target intervertebral space. The device 14 can then be pushed (e.g., by a plunger)
- out of the steering tube and into the intervertebral space. The steering tube 28 does not have to,
- 28 but can, enter the intervertebral space.
- 29 [0060] Figures 10c through 10d illustrate that the device 14 can have one or more exterior
- 30 steering rails, tracks or wires 30a and 30b, such as guidewires. The rails 30a and 30b can
- 31 slidably or fixedly and releasably engage the external surface of the segments 22 of the device

1 14. For example, the rails can pass through slots, guides, collars, cuffs or combinations thereof

- on the exterior of the segments 22. The slots, guides, collars, cuffs or combinations thereof,
- and/or the rails 30a and 30b can be coated or covered with a low-friction (e.g., PTFE) or high-
- 4 friction (e.g., knurled or toothed surface texturing) material or surface treatment or texture,
- 5 including any of the materials listed herein. The steering rails 30a and 30b can be steered or
- 6 manipulated by applying a tensile force to tensile cables within the rails, as shown by the arrows
- 7 in Figures 10d and 10e, and the flexing from Figure 10c to 10d. The rails 30a and 30b can be
- 8 pre-formed to a specific shape and can be substituted for other rails 30a and 30b that can be pre-
- 9 formed to a different shape to change the direction of delivery.
- 10 [0061] Figures 11a through 11c illustrates that the device 14 can have one or more interior
- steering rails, guide, tracks or wires 30, such as guidewires. The rails 30 can be positioned
- through the center or interior of one or more segments 22 of the device 14. The rail 20 can
- slidably or fixedly and releasably engage an internal surface, such as through a longitudinal
- 14 guide port or channel 32, of the segments 22 of the device 14. For example, ports or channels
- can extend longitudinally through the segments 22 of the device 14. The channels, and/or the
- rail 30 can be coated, covered or collared, such as with a low-friction (e.g., PTFE) or high-
- 17 friction (e.g., knurled or toothed surface texturing) material or surface treatement or texture,
- including any of the materials listed herein. The steering rail 30 can be steered or manipulated
- by applying a tensile force to tensile cables within the rail 30, as shown by the flexing from
- 20 Figure 11a to 11c. The rail 30 can be pre-formed to a specific shape and can be substituted for
- 21 one or more other rails 30 that can be pre-formed to a different shape to change the direction of
- 22 delivery.
- 23 [0062] The distal ends of the internal and/or external steering rail or rails 30 can be positioned at
- 24 the target deployment site. For example, the steering rails 30 can be placed in the intervertebral
- 25 space and can remain in the intervertebral space post-surgery, or the steering rails 30 can be
- 26 removed from the intervertebral space and the device 14 can be deployed from the rails 30 and
- 27 the device 14 can be left in the intervertebral space.
- 28 [0063] Also for example, the distal end of the steering rails 30 can be positioned at the entrance
- 29 to the intervertebral space and/or rested on the inferior and/or superior vertebral body end plate
- 30 adjacent to the target intervertebral space. The device 14 can then be pushed (e.g., by a plunger)

out of the steering rails 30 and into the intervertebral space. The steering rails 30 do not have to,

- 2 but can, enter the intervertebral space.
- 3 [0064] Figures 12a through 12f illustrate cross-sections of various rails 30, or at various lengths
- 4 along the same rail 30. Figure 12a illustrates that the cross-section of the steering rail 30 can be
- 5 circular. Figure 12b illustrates that the cross-section of the steering rail 30 can be oval. Figure
- 6 12c illustrates that the cross-section of the steering rail 30 can be multi-ovular (i.e., having a
- 7 union of two or more ovals with the same major axis). Figure 12d illustrates that the cross-
- 8 section of the steering rail 30 can be the union of rectangles intersecting at right (or another)
- 9 angle, such as a plus-sign. Figure 12e illustrates that the cross-section of the steering rail 30 can
- 10 be hexagonal. Figure 12f illustrates that the cross-section of the steering rail 30 can be
- 11 rectangular or square with sharp or rounded (chamfered) edges. The cross-section of the steering
- rail 30 can be triangular, pentagonal, heptagonal, or octagonal. The steering rail 30, whether
- internal or external to the device 14, can deliver torque around the longitudinal and/or transverse
- 14 axes of the device. The steering rail 30 can have various cross sections at various lengths along
- the rail 30. The steering rail 30 can guide, pitch, yaw and roll the device 14 into a desired
- orientation or indication. The device 14 can be delivered with one or more internal and/or
- 17 external rails 30 and/or a sheath 28 or neither.
- 18 [0065] Figure 13 illustrates a device 14 that can be attached to a deployment tool having a
- 19 controller handle 34 controllably attached to the internal steering rail 30. The internal steering
- 20 rail 30 can pass through the device 14. The steering rail 30 can be fixedly attached to the device
- 21 14 during the delivery and articulation of the device 14. The device can be steered along or
- tracking the medial surface of the Ilium 2. The device 14 can then be positioned adjacent to the
- 23 target site (e.g., the L5-S1 intervertebral space). The deployment tool can then release the device
- 24 14 from the steering rail 30 and push the device 14 into the target site.
- 25 [0066] Figures 14a and 14b illustrate that the device 14 can be delivered by being pushed along a
- steering horn, boot, or slide 36. The slide 36 can be similar to the steering tube 28, except that at
- least one wall of the slide 36 can be missing or open (e.g., the top wall is not present in the
- variation of the slide shown) compared with the steering tube 28. The missing wall can be
- completely open or replaced by one or more steering rails 30. The slide 36 can be used similar to
- 30 the steering rails 30 and/or steering tube 28. The slide 36 can be steered, flexed or articulated by

applying a tensile force to tensile cables within the rails, as shown by the arrow in Figure 14b,

- 2 and the flexing from Figure 14a to 14b.
- 3 [0067] Figures 15a and 15b illustrate that the device 14 can have six segments 22a through 22f
- 4 and five hinges 24a through 24e. The segments 22 can be attached to adjacent segments 22 by
- 5 one or more hinges, tension or steering rails or wires, screws, pins, or combinations thereof.
- 6 The hinges 24 can be pins. The segments 22 can be chained together. The segments 22 can be
- 7 identical to each other except for the distal-most segment 22a and the proximal-most segment
- 8 22f. The segments 22 or links can be box-shaped. The hinges 24, such as the pins, can be
- 9 parallel to all or some of the other hinges 24.
- 10 [0068] Figure 16 illustrates that the hinges 24 can be at acute angles to all or some of the hinges
- 11 24. The hinges 24 can be at hinge angles 38 with respect to each other. The hinge angle 38 can
- be measured between the hinge longitudinal axis 40 and the device longitudinal axis 42. The
- hinge angles 38 can be from about 80° to about 150°, more narrowly from about 90° to about
- 14 135°, yet more narrowly from about 95° to about 110°.
- 15 [0069] The device 14 can be translated and/or rotated by a handle 34 that can be removably
- attached to the device 14. The handle 34 can be screwed and/or snapped directly into the
- proximal end of the device 14, such as into the proximal-most segment 22. The handle 34 can
- compress, such as by grabbing or pinching, the proximal end of the device 14. The handle 34
- 19 can be a pusher, plunger, ram, or combinations thereof. The handle 34 and/or remainder of the
- 20 deployment tool can be rigid and/or flexible or articulatable. For example, hinged similar to the
- 21 device 14.
- 22 [0070] The segments 22 are not necessarily connected to each other by hinges. The segments 22
- can be delivered to the target site individually, or as an unattached line of segments 22.
- 24 [0071] The device 14 can be cylindrical and flexible. The implantable device 14 can be fully
- 25 flexible all the time. The device 14 can be mechanically stabilized by the deployment tool,
- steering wires, sheaths, tubes and guides. For example, the tools, wires, sheaths, tubes and
- 27 guides can provide column stability to press the device 14 into the target site (e.g., intervertebral
- 28 disc space).
- 29 [0072] The device 14 can flexible, and then locked with a tension or steering wire to stop
- 30 rotational motion of the hinges once the device is delivered to and oriented within the target site.

1 The tension wire could be tightened, for example by being tensioned by a nut to create higher

- 2 friction in each hinge 24.
- 3 [0073] Figures 17a through 17f illustrate that the device 14 can have a living hinge 44. The
- 4 living hinge 44 is a length of decreased rigidity and increased flexing within the body of the
- 5 device 14. The living hinge 44 can be formed around slots and continuous segments of
- 6 otherwise tough, durable material. The living hinge 44 can be defined be narrowing or thinning
- 7 in the body of the device 14, such that the narrowing is sufficient to provide flexibility under
- 8 reasonable torque. For example, the thickness of the unitary body of the device 14 at the living
- 9 hinge 44 can be narrowed by more than about 85%, or more than about 90%, or more than about
- 10 95%, or more than about 97%, or more than about 98.5%. The living hinge 44 can have one or
- more repeated thinnings along the length of the device 14, as shown in Figures 17a through 17f.
- 12 [0074] Figures 17a and 17b illustrate that the device 14 bends at the living hinge 44. The living
- hinges 44 can be made to control the bend and direction of the device 14. The outer surface of
- 14 the device 14 along the living hinge 44 can be smooth, for example providing low-friction
- 15 surface for sliding over bone.
- 16 [0075] Figures 17a and 17b illustrate that the living hinge 44 can be along the bottom of the
- implant device 14. Figure 17c illustrates that the living hinge 44 can be along the top of the
- device 14. Figure 17d illustrates that the living hinge 44 can be through the middle or central
- 19 axis of the device 14. Figure 17e illustrates that the living hinge 44 is discontinuous and on
- 20 opposite sides of the center of the device 44. Figure 17f illustrates that the living hinge 44 is at
- 21 an angle with respect to the longitudinal axis of the device 14, starting near the bottom of the
- device 14 and ending near the top of the device 14.
- 23 [0076] Figure 18 illustrates that the living hinge 42 can be at a non-zero angle to the central
- longitudinal axis 42 of the device 14. A first length of the living hinge 42 can be at a non-zero
- angle to a second length of the living hinge 44.
- 26 [0077] Figure 19 illustrates that the living hinge 44 can be curved. The living hinge 44 can
- 27 curve around the central longitudinal axis 42 of the device 14.
- 28 [0078] Figures 20a through 20c illustrate that the device can have three segments 22a, 22b and
- 29 22c connected by two hinges 24a and 24b. The device longitudinal axis 42 can be straight or can
- 30 have a longitudinal radius of curvature 46. The longitudinal radius of curvature 46 can be from

about 3 cm to about 100 cm, more narrowly from about 5 cm to about 20 cm, yet more narrowly

- 2 from about 7 cm to about 15 cm, for example about 15 cm, also for example about 10 cm.
- 3 [0079] The support device 14 can have a support device width 11a, a support device length 11b
- 4 and a support device height 11c. The support device width 11a can be from about 10 mm to 30
- 5 mm, or more narrowly 16 mm to about 18 mm. The support device length 11b can be from
- 6 about 30 mm to 60 mm, or more narrowly from 45 mm to about 55 mm. The support device
- 7 height 11c can be from about 1 mm to 30 mm, or more narrowly from 8 mm to about 16 mm.
- 8 [0080] The device 14 can have an anterior taper or lordosis angle 48. The taper angle 48 can be
- 9 measured between the plane of the top surface and the plane of the bottom surface of the device
- 10 14. The taper angle 48 can be from about 0° (i.e., parallel top and bottom planes) to about 45°,
- more narrowly from about 2° to about 20°, yet more narrowly from about 0° to about 12°, yet
- more narrowly from about 4° to about 10°, yet more narrowly from about 4° to about 8°, for
- example from about  $0^{\circ}$ , also for example about  $6^{\circ}$ .
- 14 [0081] The first, second, and third links or segments 22a, 22b and 22c of the flexible implantable
- device 14 may be separate or connected. One or more of the segments 22 can be rigid and/or
- 16 flexible. One or more of the segments 22 can have through-ports or segment ports 50, such as
- 17 first, second and third segment ports 50a, 50b and 50c, through the first, second and third
- segments 22a, 22b, and 22c, respectively. The segment ports 50 can extend through part of all of
- 19 the height of the respective segment 22 or the device 14 from the top to the bottom surface. One
- 20 or more of the segment ports 50 can be partially or completely filled with a bone ingrowth
- 21 matrix, bone morphogenic protein, therapeutic agents, any agent or material disclosed herein, or
- 22 combinations thereof, for example for analgesic effect or to promote bone ingrowth.
- 23 [0082] The device 14 can have a surface coating or texturing on the top, and/or bottom, and/or
- side surfaces, such as lateral teeth 52, longitudinal or angled teeth, knurling, a coating or matrix
- 25 to promote bone ingrowth, or combinations thereof.
- 26 [0083] The device 14 can have hinge teeth or knuckles 54. The hinge teeth 54 can slide by
- adjacent hinge teeth 54 to increase lateral stability during articulation and increase range of
- 28 motion (e.g., a hinge tooth 54 on one segment 22 can slide into the gap between hinge teeth 54
- 29 on the adjacent segment 22 during articulation of the device 14).

1 [0084] One or more tension and/or steering wires can be inserted and/or tensioned through guide

- 2 ports or channels 32a and 32b. The guide channels 32a and 32b can extend longitudinally
- 3 through some or all of the segments 22.
- 4 [0085] The first segment 22a and the third segment 22c can have central vertical holes 82a and
- 5 82b, respectively. The central vertical holes 82 can be attached to a deployment device, screwed
- 6 to the adjacent tissue (i.e., bone) after delivery, filled with a radiopaque material for visualization
- 7 or therapeutic or other material listed herein, or combinations thereof.
- 8 [0086] Figures 21a through 21c illustrate that device 14 can articulate. The segments 22 can
- 9 rotate with respect to each other about the hinges 24, as shown by arrows.
- 10 [0087] The first segment 22a can have a first segment longitudinal axis 42a. The second
- segment 22b can have a second segment longitudinal axis 42b. The third segment 22c can have
- 12 a third segment longitudinal axis 42c. The respective longitudinal axes can intersect at the
- adjoining hinge pins 24. The first segment longitudinal axis 42a can form a first articulation
- angle 80a with the second segment longitudinal axis 42b. The second segment longitudinal axis
- 15 42b can form a second articulation angle 80b with the third segment longitudinal axis 42c. The
- 16 first and second articulation angles 80a and 80b can be the same or different. When the device is
- in an articulated configuration, the first and/or second articulation angles 80a and/or 80b can be
- 18 from about 0° to about 90°, more narrowly from about 3° to about 75°, yet more narrowly from
- about 5° to about 60°, yet more narrowly from about 15° to about 45°.
- 20 [0088] Figures 22a through 22c illustrate that some or all of the distal-most segments 22a
- 21 through 22d can be identical. Segments 22 can be added or removed from the device 14, before
- during or after deployment to the target site, to increase or decrease the length of the device 14 to
- best fit the target site. The false hinge 24' can be a hinge component that is not attached to the
- other half of the hinge 24. The hinges 24 can snap together and apart. The articulation of each
- 25 segment 22 can be limited by the interference fit of a rotational stop 58 on the top and bottom of
- the adjacent segment 22.
- 27 [0089] The device 14 can have a deployment tool interface, such as the lateral hole 56, for
- 28 attaching to the deployment tool.
- 29 [0090] Figures 23a through 23c illustrate that a tensioning or steering wire or rail 30 can be
- 30 deployed through the channels 32 on each segment. The wire 30 can then be tensioned to
- 31 articulate and/or lock the device 14 in an articulated configuration.

1 [0091] PCT Application No. PCT/US 11/00974 filed 27 May 2011 which claims priority to U.S.

- 2 Provisional App. No. 61/349,151 filed 27 May 2010 are both herein incorporated by reference in
- 3 their entireties.
- 4 [0092] Any or all elements of the device and/or other devices or apparatuses described herein
- 5 can be made from, for example, a single or multiple stainless steel alloys, nickel titanium alloys
- 6 (e.g., Nitinol), cobalt-chrome alloys (e.g., ELGILOY® from Elgin Specialty Metals, Elgin, IL;
- 7 CONICHROME® from Carpenter Metals Corp., Wyomissing, PA), nickel-cobalt alloys (e.g.,
- 8 MP35N® from Magellan Industrial Trading Company, Inc., Westport, CT), molybdenum alloys
- 9 (e.g., molybdenum TZM alloy, for example as disclosed in International Pub. No. WO
- 10 03/082363 A2, published 9 October 2003, which is herein incorporated by reference in its
- entirety), tungsten-rhenium alloys, for example, as disclosed in International Pub. No. WO
- 12 03/082363, polymers such as polyethylene teraphathalate (PET)/polyester (e.g., DACRON®
- from E. I. Du Pont de Nemours and Company, Wilmington, DE), polypropylene, (PET),
- 14 polytetrafluoroethylene (PTFE), expanded PTFE (ePTFE), polyether ketone (PEK), polyether
- ether ketone (PEEK), poly ether ketone ketone (PEKK) (also poly aryl ether ketone ketone),
- 16 nylon, polyether-block co-polyamide polymers (e.g., PEBAX® from ATOFINA, Paris, France),
- 17 aliphatic polyether polyurethanes (e.g., TECOFLEX® from Thermedics Polymer Products,
- Wilmington, MA), polyvinyl chloride (PVC), polyurethane, thermoplastic, fluorinated ethylene
- 19 propylene (FEP), absorbable or resorbable polymers such as polyglycolic acid (PGA), polylactic
- 20 acid (PLA), polycaprolactone (PCL), polyethyl acrylate (PEA), polydioxanone (PDS), and
- 21 pseudo-polyamino tyrosine-based acids, extruded collagen, silicone, zinc, echogenic, radioactive,
- radiopaque materials, a biomaterial (e.g., cadaver tissue, collagen, allograft, autograft, xenograft,
- 23 bone cement, morselized bone, osteogenic powder, beads of bone) any of the other materials
- 24 listed herein or combinations thereof. Examples of radiopaque materials are barium sulfate, zinc
- 25 oxide, titanium, stainless steel, nickel-titanium alloys, tantalum and gold.
- 26 [0093] Any or all elements of the device and/or other devices or apparatuses described herein.
- 27 can be, have, and/or be completely or partially coated with agents and/or a matrix a matrix for
- cell ingrowth or used with a fabric, for example a covering (not shown) that acts as a matrix for
- 29 cell ingrowth. The matrix and/or fabric can be, for example, polyester (e.g., DACRON® from E.
- 30 I. Du Pont de Nemours and Company, Wilmington, DE), polypropylene, PTFE, ePTFE, nylon,
- 31 extruded collagen, silicone or combinations thereof.

1 [0094] The device and/or elements of the device and/or other devices or apparatuses described

- 2 herein and/or the fabric can be filled, coated, layered and/or otherwise made with and/or from
- 3 cements, fillers, glues, and/or an agent delivery matrix known to one having ordinary skill in the
- 4 art and/or a therapeutic and/or diagnostic agent. Any of these cements and/or fillers and/or glues
- 5 can be osteogenic and osteoinductive growth factors.
- 6 [0095] Examples of such cements and/or fillers includes bone chips, demineralized bone matrix
- 7 (DBM), calcium sulfate, coralline hydroxyapatite, biocoral, tricalcium phosphate, calcium
- 8 phosphate, polymethyl methacrylate (PMMA), biodegradable ceramics, bioactive glasses,
- 9 hyaluronic acid, lactoferrin, bone morphogenic proteins (BMPs) such as recombinant human
- bone morphogenetic proteins (rhBMPs), other materials described herein, or combinations
- 11 thereof.
- 12 [0096] The agents within these matrices can include any agent disclosed herein or combinations
- 13 thereof, including radioactive materials; radiopaque materials; cytogenic agents; cytotoxic
- 14 agents; cytostatic agents; thrombogenic agents, for example polyurethane, cellulose acetate
- 15 polymer mixed with bismuth trioxide, and ethylene vinyl alcohol; lubricious, hydrophilic
- materials; phosphor cholene; anti-inflammatory agents, for example non-steroidal anti-
- inflammatories (NSAIDs) such as cyclooxygenase-1 (COX-1) inhibitors (e.g., acetylsalicylic
- acid, for example ASPIRIN® from Bayer AG, Leverkusen, Germany; ibuprofen, for example
- 19 ADVIL® from Wyeth, Collegeville, PA; indomethacin; mefenamic acid), COX-2 inhibitors
- 20 (e.g., VIOXX® from Merck & Co., Inc., Whitehouse Station, NJ; CELEBREX® from
- 21 Pharmacia Corp., Peapack, NJ; COX-1 inhibitors); immunosuppressive agents, for example
- 22 Sirolimus (RAPAMUNE®, from Wyeth, Collegeville, PA), or matrix metalloproteinase (MMP)
- 23 inhibitors (e.g., tetracycline and tetracycline derivatives) that act early within the pathways of an
- 24 inflammatory response. Examples of other agents are provided in Walton et al, Inhibition of
- 25 Prostoglandin E<sub>2</sub> Synthesis in Abdominal Aortic Aneurysms, Circulation, July 6, 1999, 48-54;
- 26 Tambiah et al. Provocation of Experimental Aortic Inflammation Mediators and Chlamydia
- 27 Pneumoniae, Brit. J. Surgery 88 (7), 935-940; Franklin et al, Uptake of Tetracycline by Aortic
- Aneurysm Wall and Its Effect on Inflammation and Proteolysis, Brit. J. Surgery 86 (6), 771-775;
- 29 Xu et al, Sp1 Increases Expression of Cyclooxygenase-2 in Hypoxic Vascular Endothelium, J.
- 30 Biological Chemistry 275 (32) 24583-24589; and Pyo et al, Targeted Gene Disruption of Matrix
- 31 Metalloproteinase-9 (Gelatinase B) Suppresses Development of Experimental Abdominal Aortic

1 Aneurysms, J. Clinical Investigation 105 (11), 1641-1649 which are all incorporated by

2 reference in their entireties.

3

### 4 **[0097]** METHODS OF USE

- 5 [0098] Figure 24 illustrates that a straight or curved transosseous delivery channel 94 can be
- 6 drilled, chiseled, punched, or a combination thereof, through the iliac bone 2 and/or the sacral ala
- 7 90, for example passing through the sacroiliac joint 92. The transosseous delivery channel 94
- 8 have a first length or first channel through the iliac 2 and a second length or second channel
- 9 through the sacrum S1. The first length of the transosseous delivery channel 94 can be aligned
- with the second length of the transosseous delivery channel 94, for example to form a
- 11 substantially continuous channel. The transosseous delivery channel 94 can have a laterally-
- located channel entry port 96 laterally outside of the sacral ala 90 and/or iliac bone 2. The
- transosseous delivery channel 94 can have a channel exit port 98 adjacent to the L5-S1
- intervertebral disc space. For example, the channel exit port 98 can be in the S1 endplate. The
- channel exit port 98 can be positioned so the circumference of the channel exit port 98
- tangentially coincides with or is closely adjacent to (e.g., within about 2 cm, more narrowly
- within about 1 cm, more narrowly within about 5 mm, yet more narrowly within about 2 mm)
- with the edge of the S1 vertebral endplate 100.
- 19 [0100] The L5-S1 intervertebral space can be partially or completely void of soft tissue, as
- shown, for example from a discectomy performed before insertion of the support device 14. For
- 21 example, the discectomy can be performed by the method and device shown in U.S. Provisional
- 22 Patent Application No. 61/526,630 filed 23 August 2011, which is incorporated by reference
- 23 herein in its entirety.
- 24 [0101] Figure 25 illustrates that the support device 14 can be inserted, as shown by arrow
- 25 102, medially through the channel entry port 96 of the transosseous delivery channel 94. The
- device 14 can be removably and/or articulatably attached to a deployment tool shaft 104.
- 27 [0102] Figure 26 illustrates that the shaft 104 can be further translated, as shown by arrow
- 28 106, into the transosseous delivery channel 94. The support device 14 can translate toward and
- 29 into the L5-S1 intervertebral disc space. The distal tip of the support device 14 can enter the L5-
- 30 S1 intervertebral disc. The support device 14 can enter the target site of the L5-S1 intervertebral

disc directly from the transosseous delivery channel 94 without passing through any soft tissue


- 2 between the L5-S1 intervertebral disc and the iliac bone 2.
- 3 [0103] Figure 27 illustrates that the shaft 104 can be further translated, as shown by arrow
- 4 108, medially through the transosseous delivery channel 94. The device 14 can translate, as
- 5 shown by arrow 110, through the L5-S1 intervertebral disc space and the L5 and/or the S1
- 6 vertebra. The support device 14 can articulate, as shown by arrow 112. One or more of the
- 7 hinges 24 can rotate, articulating the segments 22. The hinges 24 can be controllably rotatably
- 8 locked and unlocked, for example, by controls on the handle of the deployment tool (of which
- 9 the shaft 104 is a part).
- 10 [0104] The support device 14 can then be further translated, such as being pushed and/or
- vibrated (e.g., manually, ultrasonically), for example, medially and laterally, and/or superior and
- inferiorly, and/or anteriorly and posteriorly. The through ports and/or cavities and/or recesses 50
- in the device 14 can be partially and/or completely filled bone morphogenic protein, therapeutic
- 14 agents, other materials listed herein or combinations thereof. The support device 14 can deliver
- a cauterizing electrical energy from the deployment tool. The support device 14 and shaft 104
- can have one or more longitudinal lumens that can be used to irrigate (e.g., with analgesic agents,
- 17 saline, anesthetic agents, bone morphogenic proteins, visualization agents, other agents described
- herein, or combinations thereof) and/or aspirate (e.g., to remove irrigated material and/or debris)
- 19 the target site (e.g., the L5-S1 intervertebral disc space).
- 20 [0105] Figure 28 illustrates that before, during or after the support device 14 is positioned in
- 21 the L5-S1 intervertebral space, the shaft 104 can detach from the support device 14 and be
- translated laterally, as shown by arrow 114, from the L5-S1 intervertebral disc space and the
- 23 transosseous delivery channel 94. The deployment tool shaft 104 can remove or reposition the
- support device 14, or leave the support device 14 in place in the L5-S1 space.
- 25 [0106] The method shown in Figures 25 through 28 can be repeated to deliver multiple
- support devices 14 to one or more intervertebral spaces. For example, one, two, three or more
- 27 support devices 14 can be positioned in the L4-L5 intervertebral space and/or the L5-S1
- 28 intervertebral space. The one, two, three or more support devices 14 positioned in the L4-L5
- 29 and/or L5-S1 intervertebral spaces, can mechanically support the adjacent vertebrae and/or fix
- 30 the adjacent vertebrae to each other. Bone ingrowth can occer through the through ports 50, for
- 31 example fusing the support devices 14 to the respective surrounding vertebrae.


1 [0107] Figures 29 through 31 illustrate that the transosseous delivery channel 94 can have a

- 2 coronal delivery angle 120 measured to the coronal plane 122, a sagittal delivery angle 124
- 3 measured to the sagittal plane 126, and a transverse delivery angle 128 measured to the
- 4 transverse plane 130. The coronal delivery angle 120 can be from about 0° to about 25°, for
- 5 example about 12°. The sagittal delivery angle 124 can be from about 65° to about 90°, for
- example about 75°. The transverse delivery angle 128 can be from about  $0^{\circ}$  to about  $20^{\circ}$ , for
- 7 example about 10°. The support device 14 and shaft 104 can be configured so the support device
- 8 14 can exit the channel exit port 98 (e.g., directly into the L5-S1 intervertebral disc) and
- 9 articulate sufficiently to enter and pass through all or a significant portion (e.g., more than about
- 10 40%, yet more narrowly more than about 50%, yet more narrowly more than about 75%) of the
- 11 width of the L5-S1 intervertebral space.
- 12 [0108] Figures 30 and 31 show one or both femurs 132 for illustrative purposes.
- 13 [0109] Figures 32a through 32d, and separately Figures 33a through 33c illustrate the
- deployment of the support device 14 into the L5-S1 intervertebral disc space target site, as
- described for Figures 24-27. The support device 14 can be delivered to a complete or partial
- discectomy target site 138 in the L5-S1 space.
- 17 [0110] Figure 33d illustrates the shaft 104 can be rotated, as shown by arrow 134, about the
- 18 longitudinal axis of the shaft before during or after the support device 14 is positioned in the L5-
- 19 S1 intervertebral disc space target site. The support device 14 can rotate, as shown by arrow
- 20 136, in the L5-S1 intervertebral disc space, for example to control and position the support
- 21 device 14 to an angular orientation in the transverse plane 130.
- 22 [0111] Any elements described herein as singular can be pluralized (i.e., anything described
- as "one" can be more than one). Any species element of a genus element can have the
- 24 characteristics or elements of any other species element of that genus. The above-described
- 25 configurations, elements or complete assemblies and methods and their elements for carrying out
- 26 the invention, and variations of aspects of the invention can be combined and modified with each
- 27 other in any combination.

1	CLAIMS
2	We claim:
3	1. A method for inserting an implant to a target site between a first vertebra and a second
4	vertebra comprising:
5	creating a first channel through the ilium;
6	creating a second channel through the sacrum, wherein the first channel is aligned with
7	the second channel;
8	inserting a first rigid section of the implant through the first channel and the second
9	channel into the target site,
10	rotating a second rigid section of the implant with respect to the first rigid section,
11	wherein the first rigid section is hingedly attached to the second rigid section; and
12	inserting the second rigid section of the implant into the target site.
13	
14	2. The method of Claim 1, wherein the second channel passes through a vertebral endplate.
15	
16	3. A biological implant support device for providing orthopedic support comprising:
17	a first rigid section at a first terminal end of the device, the first rigid section having a
18	first top plate and a first bottom plate; and
19	a second rigid section having a second top plate and a second bottom plate;
20	wherein the first rigid section is rotatably attached to the second rigid section, and
21	wherein the top and bottom plates are configured to interface with hard tissue.
22	
23	4. A method for inserting a support device to a target site in a spine adjacent to a first vertebra
24	comprising;
25	creating a channel through a non-vertebral bone and a vertebral end plate;
26	inserting a first rigid section of the device through the channel and into the target site,
27	inserting a second rigid section of the device through the channel,
28	rotating the second rigid section of the implant with respect to the first rigid section,
29	wherein the first rigid section is hingedly attached to the second rigid section; and
30	inserting the second rigid section of the implant into the target site.


5. The method of Claim 4, wherein creating the channel comprises drilling with a flexible drill. 1 2 6. The method of Claim 4, wherein the non-vertebral bone comprises the pelvis. 3 4 7. The method of Claim 4, wherein the non-vertebral bone comprises the ilium. 5 6 8. The method of Claim 4, wherein the non-vertebral bone comprises the sacrum. 7 8 9. A method for inserting a support device to a target site in a spine adjacent a first vertebra 9 comprising: 10 creating a channel through a non-vertebral bone and a vertebra; and 11 inserting the support device through the channel and into the target site. 12 13 10. A biological implant device for providing orthopedic support comprising: 14 a first rigid section with a first top plate and a first bottom plate; and 15 a second rigid section with a second top plate and a second bottom plate; 16 17 wherein the first rigid section is rotatably attached to the second rigid section, and wherein the top and bottom plates are configured to interface with hard tissue. 18







## **NOT INVENTION**


Fig. 3





NOT INVENTION Fig. 4














Fig. 11a




Fig. 11b





Fig. 11c



30



Fig. 12a

Fig. 12b

Fig. 12c







Fig. 12d

Fig. 12e

Fig. 12f






Fig. 13










Fig. 15b








Fig. 17b




Fig. 17c




Fig. 17d




Fig. 17e




Fig. 17f




Fig. 18









Fig. 19























Fig. 24




Fig. 25




Fig. 26




Fig. 27




Fig. 28





Fig. 29






Fig. 31

