

US 20060205040A1

# (19) United States (12) Patent Application Publication (10) Pub. No.: US 2006/0205040 A1

# (10) Pub. No.: US 2006/0205040 A1 (43) Pub. Date: Sep. 14, 2006

# Sampath

(54) COMPOSITIONS FOR USE IN IDENTIFICATION OF ADVENTITIOUS VIRUSES

(76) Inventor: **Rangarajan Sampath**, San Diego, CA (US)

Correspondence Address: MEDLEN & CARROLL LLP 101 HOWARD STREET SUITE 350 SAN FRANCISCO, CA 94105 (US)

- (21) Appl. No.: **11/368,233**
- (22) Filed: Mar. 3, 2006

**Related U.S. Application Data** 

(60) Provisional application No. 60/658,248, filed on Mar.
3, 2005. Provisional application No. 60/705,631, filed on Aug. 3, 2005. Provisional application No. 60/732, 539, filed on Nov. 1, 2005.

#### **Publication Classification**

- (51) Int. Cl. *C12Q* 1/68 (2006.01) *C12P* 19/34 (2006.01)
- (52) U.S. Cl. ..... 435/91.1; 435/6

# (57) **ABSTRACT**

The present invention provides compositions, kits and methods for rapid identification and quantification of adventitious contaminant viruses by molecular mass and base composition analysis.





Figure 1





# RELATED APPLICATIONS

[0001] The present application 1) claims the benefit of priority to U.S. Provisional Application Ser. No. 60/658,248, filed Mar. 3, 2005; 2) claims the benefit of priority to U.S. Provisional Application Ser. No. 60/705,631, filed Aug. 3, 2005; 3) claims the benefit of priority to U.S. Provisional Application Ser. No. 60/732,539, filed Nov. 1, 2005 and 4) claims the benefit of priority to U.S. Provisional Application Ser. No. 60/740,617, filed Nov. 28, 2005. Each of the above listed U.S. Provisional Applications is incorporated herein by reference in entirety. Methods disclosed in U.S. application Ser. Nos. 10/156,608, 09/891,793, 10/418,514, 10/660, 997, 10/660,122, 10,660,996, 10/660,998, 10/728,486, 10/405,756, 11/060,135, and 11/073,362, are commonly owned and incorporated herein by reference in their entirety for any purpose.

#### STATEMENT OF GOVERNMENT SUPPORT

**[0002]** This invention was made with United States Government support under NIH contract N01 Al40100. The United States Government may have certain rights in the invention.

#### FIELD OF THE INVENTION

**[0003]** The present invention provides compositions, kits and methods for rapid identification and quantification of adventitious contaminant viruses by molecular mass and base composition analysis.

#### BACKGROUND OF THE INVENTION

#### A. Adventitious Viruses

[0004] Adventitious viruses represent a major risk associated with the use of cell-substrate derived biologicals, including vaccines and antibodies, for human use. The possibility for viral contamination exists in primary cultures and established cultures, as well as Master Cell Banks, end-of-production cells, and bulk harvest fluids. This is a major obstacle to the use of neoplastic-immortalized cells for which the mechanism of transformation is unknown is that these could have a higher risk of containing oncogenic viruses. Extensive testing for the presence of potential extraneous agents is therefore required to ensure the safety of the vaccines. Among the methods used for this purpose are animal inoculations, electron microscopy and in vitro molecular and antibody assays that provide a screen for viral agents. Another critical consideration for assessing the safety concerns associated with viral vaccines is the detection of endogenous retroviral sequences while using avian, murine, non-human primate, and human cell lines. Endogenous retroviral sequences are an integral part of eukaryotic genomes, and while the majority of these sequences are defective, a few can produce infectious virus, either spontaneously upon long-term culture. These can also be induced upon treatment with various chemical or other agents that may be part of the normal production system. The activation of an endogenous, infectious retrovirus in a cell substrate that is used for the production of biologics is an important safety concern, especially in the case of live, viral vaccines, where minimal purification and inactivation steps are used in order to preserve high vaccine potency.

[0005] The currently established methods for measuring RT-activity include the highly sensitive, product-enhanced reverse transcriptase assays (PERT) that can detect 1-10 virions and transmission electron microscopy (TEM) to analyze infective retroviruses particles. However, the above techniques are not specific and do not provide any information regarding the source of the RT activity. PCR-based detection of retroviruses can be used in combination with other assays such as reverse transcriptase, electron microscopy infectivity or co-cultivation to increase the sensitivity of detection or to identify a particular adventitious agent present in the test sample. Further, while some studies demonstrate that a low level of RT activity is not generally associated with a replicating agent; major concerns remain regarding the consequences of the presence of such nonproductive, non-replicating defective infections in the vaccine, as there is the potential for integration into the host genome.

[0006] Retrovirus-induced tumorigenesis can involve the generation of a novel pathogenic virus by recombination between replication-competent and -defective sequences and/or activation of a cellular oncogene by a long terminal repeat (LTR) due to upstream or downstream insertion of retrovirus sequences. To address the possible integration of extraneous retroviral sequences in human cells by RT-containing particles, multiple PCR strategies have been used. These include direct PCR of DNase-treated inoculum using primers from the highly conserved pol region and Alu PCR using LTR primers in conjunction with Alu primers that specifically amplify viral-cellular DNA junctions of integrants.

**[0007]** Future strategies to detect adventitious agents must address three fundamental problems. First, there are large numbers of known viral agents that are potential contaminants, each with a large number of potential strain variants. Second, history has shown that not all adventitious agents fall into anticipated families of viruses, so unanticipated virus families must also be considered. Third, the test must be practical to perform on a large number of samples in a standardized, high-throughput, quality-controlled fashion. The premise of this proposal is that we can leverage recently developed and validated methods using mass spectrometry analysis of broad-range PCR reactions for rapid, sensitive, cost-effective detection of broad ranges of adventitious agents, including previously unknown/uncharacterized viruses and endogenous retroviruses.

### B. Drug Resistance

**[0008]** Drug resistance in bacteria and viruses is frequently mediated by point mutations in key genes whose gene products interact directly or indirectly with the drug. While there are several methods available for identification of single nucleotide polymorphisms (SNPs) in nucleic acid sequences, the functional unit that encodes each amino acid is the codon, where three successive nucleotides are responsible for encoding each amino acid. Mutations in any of the three nucleotides may or may not result in a mutation in the encoded amino acid, depending upon the particular amino acid and the rules of the genetic code. Because the genetic code is deciphered as a sequence, both the identity and the order of the nucleotides are important in determining the encoded amino acid. Thus, DNA sequencing has become the method of choice for analysis of mutations that result in amino acid changes. DNA sequencing has significant disadvantages as an analysis method for routine use a clinical laboratory setting. It is still relatively expensive and labor intensive, and thus is used only for very important analyses. An example of this is determination of drug resistance in viruses such as HIV and in bacteria such as methicillinresistant Staphylococcus aureus (MRSA). Drug resistance in HIV has now emerged as a significant problem in both untreated and drug-treated patient populations. The decision to select a particular drug-treatment regimen that the virus will respond to is critical to success of therapy. Drug resistance testing has been shown to improve the clinical outcome in HIV-infected individuals and thus is now recommended for new infections or for patients infected as long as two years or more prior to initiating therapy, in the case of antiretroviral failures and during pregnancy. Thus, despite the costs, DNA sequencing is currently being used for determination of viral drug resistance. Typically, a serum sample is analyzed by PCR amplification of the reverse transcriptase and protease genes, followed by sequencing of approximately 900 nucleotides of the reverse transcriptase gene and 300 nucleotides of the protease gene. The DNA sequence is then used to determine the optimal drug regimen. A drawback of sequencing is that DNA sequencing technology for identification of drug-resistant viruses is that it is not easily able to identify the components present in a mixed sample, particularly in a scenario where a fraction of the virus population has mutated. DNA sequencing was developed on the assumption that the sample being analyzed is homogeneous. However, the HIV populations that infect humans are not homogeneous, and RNA viruses such as HIV are known to rapidly mutate, creating a population of mixed sequences in each infected individual. In the presence of drug selection, mutations that mediate drug resistance that occur at low frequency grow with a selective advantage and eventually can dominate the population, causing treatment failure. In this scenario, the mutant virus starts out as an undetectable fraction of the population which increases to a higher percentage over time. It would be valuable to identify drug resistant virus populations early, before they have a chance to increase the viral load. DNA sequencing methods can identify mixed populations, but do so poorly. In a recent publication using the ABI PRISM 3100 genetic analyzer, it was reported that a viral mixture containing approximately

was reported that a viral mixture containing approximately 40% of the mutant viral population can be detected with 95% confidence. However, 40% of a typical viral load (1,800 to 10,500 HIV copies/ml) means a blood burden (assuming 5 liters of blood) of up to 21 million drugresistant viral copies. Other analytical methods are capable of identifying mutations with more sensitivity than sequencing, but these methods are time consuming, laborious and not amenable to high throughput processes.

**[0009]** Thus, there is a need for rapid and cost effective methods that can be applied as alternatives to sequencing in genomic analysis for variations that mediate amino acid changes. The present invention satisfies this need. The present invention provides, inter alia, methods of identifying adventitious contaminant viruses. Also provided are oligonucleotide primers, compositions and kits containing the oligonucleotide primers, which produce amplification products whose molecular masses provide the means to identify adventitious contaminant viruses at the sub-species level.

### SUMMARY OF THE INVENTION

**[0010]** The present invention provides compositions, kits and methods for rapid identification and quantification of adventitious contaminant viruses by molecular mass and base composition analysis.

**[0011]** One embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 47.

**[0012]** Another embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 286.

**[0013]** Another embodiment is a composition of is an oligonucleotide primer pair including an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 47 and an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 286.

**[0014]** One embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 70.

**[0015]** Another embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 286.

**[0016]** Another embodiment is a composition of is an oligonucleotide primer pair including an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 70 and an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 286.

**[0017]** One embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 165.

**[0018]** Another embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 286.

**[0019]** Another embodiment is a composition of is an oligonucleotide primer pair including an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 165 and an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 286.

**[0020]** One embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 122.

**[0021]** Another embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 275.

**[0022]** Another embodiment is a composition of is an oligonucleotide primer pair including an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 122 and an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 275.

**[0023]** One embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 100.

**[0024]** Another embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 336.

**[0025]** Another embodiment is a composition of is an oligonucleotide primer pair including an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 100 and an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 336.

**[0026]** One embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 61.

**[0027]** Another embodiment is an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 324.

**[0028]** Another embodiment is a composition of is an oligonucleotide primer pair including an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 61 and an oligonucleotide primer 14 to 35 nucleobases in length having at least 70% sequence identity with SEQ ID NO: 324.

**[0029]** In some embodiments, either or both of the primers of the primer pair contain at least one modified nucleobase such as 5-propynyluracil or 5-propynylcytosine for example.

**[0030]** In some embodiments, either or both of the primers of the primer pair comprises at least one universal nucleobase such as inosine for example.

**[0031]** In some embodiments, either or both of the primers of the primer pair comprises at least one non-templated T residue on the 5'-end.

**[0032]** In some embodiments, either or both of the primers of the primer pair comprises at least one non-template tag.

**[0033]** In some embodiments, either or both of the primers of the primer pair comprises at least one molecular mass modifying tag.

**[0034]** Some embodiments are kits that contain the primer pair compositions. In some embodiments, each member of the one or more primer pairs of the kit is of a length of 14 to 35 nucleobases and has 70% to 100% sequence identity with the corresponding member from the group of primer pairs represented by SEQ ID NOs: 70:286, 165:286, 122:275, 100:336, and 61:324.

**[0035]** Some embodiments of the kits contain at least one calibration polynucleotide.

**[0036]** Some embodiments of the kits contain at least one anion exchange functional group linked to a magnetic bead.

[0037] In some embodiments, the present invention provides primers and compositions comprising pairs of primers, and kits containing the same, and methods for use in identification of adventitious contaminant viruses. The primers are designed to produce amplification products of DNA encoding genes that have conserved and variable regions across a given viral family. The invention further provides compositions comprising pairs of primers and kits containing the same, which are designed to provide species and sub-species characterization of adventitious contaminant viruses. **[0038]** In some embodiments, the present invention also provides methods for identification of adventitious contaminant viruses. Nucleic acid from the virus is amplified using the primers described above to obtain an amplification product. The molecular mass of the amplification product is measured. Optionally, the base composition of the amplification product is determined from the molecular mass. The molecular mass or base composition is compared with a plurality of molecular masses or base compositions of known adventitious contaminant virus identifying amplicons, wherein a match between the molecular mass or base composition and a member of the plurality of molecular masses or base compositions identifies the adventitious contaminant virus. In some embodiments, the molecular mass is measured by mass spectrometry.

[0039] In some embodiments, the present invention is also directed to a method for determining the presence or absence of an adventitious contaminant virus in a sample. Nucleic acid from the sample is amplified using the composition described above to obtain an amplification product. The molecular mass of the amplification product is determined. Optionally, the base composition of the amplification product is determined from the molecular mass. The molecular mass or base composition of the amplification product is compared with the known molecular masses or base compositions of one or more known adventitious contaminant virus identifying amplicons, wherein a match between the molecular mass or base composition of the amplification product and the molecular mass or base composition of one or more known adventitious contaminant virus identifying amplicons indicates the presence of the adventitious contaminant virus in the sample. In some embodiments, the molecular mass is measured by mass spectrometry.

[0040] In some embodiments, the present invention also provides methods for determination of the quantity of an unknown adventitious contaminant virus in a sample. The sample is contacted with the composition described above and a known quantity of a calibration polynucleotide comprising a calibration sequence. Nucleic acid from the unknown adventitious contaminant virus in the sample is concurrently amplified with the composition described above and nucleic acid from the calibration polynucleotide in the sample is concurrently amplified with the composition described above to obtain a first amplification product comprising an adventitious contaminant virus identifying amplicon and a second amplification product comprising a calibration amplicon. The molecular mass and abundance for the adventitious contaminant virus identifying amplicon and the calibration amplicon is determined. The adventitious contaminant virus identifying amplicon is distinguished from the calibration amplicon based on molecular mass, wherein comparison of adventitious contaminant virus identifying amplicon abundance and calibration amplicon abundance indicates the quantity of adventitious contaminant virus in the sample. In some embodiments, the base composition of the adventitious contaminant virus identifying amplicon is determined.

**[0041]** In some embodiments, the present invention provides methods for detecting or quantifying adventitious contaminant virus by combining a nucleic acid amplification process with a mass determination process. In some embodiments, such methods identify or otherwise analyze the adventitious contaminant virus by comparing mass informa-

tion from an amplification product with a calibration or control product. Such methods can be carried out in a highly multiplexed and/or parallel manner allowing for the analysis of as many as 300 samples per 24 hours on a single mass measurement platform. The accuracy of the mass determination methods in some embodiments of the present invention permits allows for the ability to discriminate between different adventitious viruses such as members of the following families: p

[0042] Papillomaviridae, Polyomaviridae, Retroviridae, Parvoviridae, Herpesviridae (Human herpesviruses 1 through 8, Bovine herpesvirus, Canine herpesvirus and Simian cytomegalovirus), Hepadnaviridae (Hepatitis B virus), Hepeviridae (Hepatitis E virus), Deltavirus (Hepatitis delta virus), Adenoviridae (Human adenoviruses A-F and murine adenovirus), Flaviviridae (Bovine viral diarrhea virus, TBE, Yellow fever virus, Dengue viruses 1-4, WNV and hepatitis C virus), Paramyxoviridae (Pneumonia virus of mice, Sendai virus, and Simian parainfluenza virus 5), Togaviridae (Western equine encephalomyelitis virus), Picornaviridae (Polio (types 1-13), Human hepatitis A, Human coxsackievirus, Human cardiovirus, Human rhinovirus and Bovine rhinovirus), Reoviridae (Mouse rotavirus, reovirus type 3 and Colorado tick fever virus), and Rhabdoviridae (vesicular stomatitis virus).

#### BRIEF DESCRIPTION OF THE DRAWINGS

**[0043]** The foregoing summary of the invention, as well as the following detailed description of the invention, is better understood when read in conjunction with the accompanying drawings which are included by way of example and not by way of limitation.

**[0044] FIG. 1**: process diagram illustrating a representative primer pair selection process.

**[0045]** FIG. 2 is a process diagram illustrating an embodiment of the calibration method.

### DEFINITIONS

**[0046]** As used herein, the term "abundance" refers to an amount. The amount may be described in terms of concentration which are common in molecular biology such as "copy number,""pfu or plate-forming unit" which are well known to those with ordinary skill. Concentration may be relative to a known standard or may be absolute.

**[0047]** As used herein an "adventitious virus" or "adventitious viral agent" refers to a virus contaminant present within a biological product, including, for example, vaccines, cell lines and other cell-derived products. In some cases, the biological product may provide a favorable environment for the survival of the virus. In some embodiments, the biological products are those useful in various experimental conditions for research in biotechnology and clinical diagnosis or treatment in pharmacology.

**[0048]** As used herein, the term "amplifiable nucleic acid" is used in reference to nucleic acids that may be amplified by any amplification method. It is contemplated that "amplifiable nucleic acid" also comprises "sample template."

**[0049]** As used herein the term "amplification" refers to a special case of nucleic acid replication involving template specificity. It is to be contrasted with non-specific template

replication (i.e., replication that is template-dependent but not dependent on a specific template). Template specificity is here distinguished from fidelity of replication (i.e., synthesis of the proper polynucleotide sequence) and nucleotide (ribo- or deoxyribo-) specificity. Template specificity is frequently described in terms of "target" specificity. Target sequences are "targets" in the sense that they are sought to be sorted out from other nucleic acid. Amplification techniques have been designed primarily for this sorting out. Template specificity is achieved in most amplification techniques by the choice of enzyme. Amplification enzymes are enzymes that, under conditions they are used, will process only specific sequences of nucleic acid in a heterogeneous mixture of nucleic acid. For example, in the case of  $Q\beta$ replicase, MDV-1 RNA is the specific template for the replicase (D. L. Kacian et al., Proc. Natl. Acad. Sci. USA 69:3038 [1972]). Other nucleic acid will not be replicated by this amplification enzyme. Similarly, in the case of T7 RNA polymerase, this amplification enzyme has a stringent specificity for its own promoters (Chamberlin et al., Nature 228:227 [1970]). In the case of T4 DNA ligase, the enzyme will not ligate the two oligonucleotides or polynucleotides, where there is a mismatch between the oligonucleotide or polynucleotide substrate and the template at the ligation junction (D. Y. Wu and R. B. Wallace, Genomics 4:560 [1989]). Finally, Taq and Pfu polymerases, by virtue of their ability to function at high temperature, are found to display high specificity for the sequences bounded and thus defined by the primers; the high temperature results in thermodynamic conditions that favor primer hybridization with the target sequences and not hybridization with non-target sequences (H. A. Erlich (ed.), PCR Technology, Stockton Press [1989]).

**[0050]** As used herein, the term "amplification reagents" refers to those reagents (deoxyribonucleotide triphosphates, buffer, etc.), needed for amplification, excluding primers, nucleic acid template, and the amplification enzyme. Typically, amplification reagents along with other reaction components are placed and contained in a reaction vessel (test tube, microwell, etc.).

**[0051]** As used herein, the term "anion exchange functional group" refers to a positively charged functional group capable of binding an anion through an electrostatic interaction. The most well known anion exchange functional groups are the amines, including primary, secondary, tertiary and quaternary amines.

**[0052]** The term "bacteria" or "bacterium" refers to any member of the groups of eubacteria and archaebacteria.

**[0053]** As used herein, a "base composition" is the exact number of each nucleobase (for example, A, T, C and G). For example, amplification of nucleic acid of *Neisseria meningitidis* with a primer pair that produces an amplification product from nucleic acid of 23S rRNA that has a molecular mass (sense strand) of 28480.75124, from which a base composition of A25 G27 C22 T18 is assigned from a list of possible base compositions calculated from the molecular mass using standard known molecular masses of each of the four nucleobases.

**[0054]** As used herein, a "base composition probability cloud" is a representation of the diversity in base composition resulting from a variation in sequence that occurs among different isolates of a given species. The "base

composition probability cloud" represents the base composition constraints for each species and is typically visualized using a pseudo four-dimensional plot.

**[0055]** In the context of this invention, a "bioagent" is any organism, cell, or virus, living or dead, or a nucleic acid derived from such an organism, cell or virus. Examples of bioagents include, but are not limited, to cells, (including but not limited to human clinical samples, bacterial cells and other pathogens), viruses, fungi, protists, parasites, and pathogenicity markers (including but not limited to: pathogenicity islands, antibiotic resistance genes, virulence factors, toxin genes and other bioregulating compounds). Samples may be alive or dead or in a vegetative state (for example, vegetative bacteria or spores) and may be encapsulated or bioengineered. In the context of this invention, a "pathogen" is a bioagent which causes a disease or disorder.

**[0056]** As used herein, a "bioagent division" is defined as group of bioagents above the species level and includes but is not limited to, orders, families, classes, clades, genera or other such groupings of bioagents above the species level.

**[0057]** As used herein, the term "bioagent identifying amplicon" refers to a polynucleotide that is amplified from a bioagent in an amplification reaction and which 1) provides sufficient variability to distinguish each individual bioagent and 2) whose molecular mass is amenable to molecular mass determination.

**[0058]** As used herein, the term "biological product" refers to any product originating from an organism. Biological products are often products of processes of biotechnology. Examples of biological products include, but are not limited to: cultured cell lines, cellular components, antibodies, proteins and other cell-derived biomolecules, growth media, growth harvest fluids, natural products and biopharmaceutical products.

**[0059]** The terms 'biowarfare agent" and "bioweapon" are synonymous and refer to a bacterium, virus, fungus or protozoan that could be deployed as a weapon to cause bodily harm to individuals by military or terrorist groups.

**[0060]** In context of this invention, the term "broad range survey primer pair" refers to a primer pair designed to produce bioagent identifying amplicons across different broad groupings of bioagents. For example, the ribosomal RNA-targeted primer pairs are broad range survey primer pairs.

**[0061]** The term "calibration amplicon" refers to a nucleic acid segment representing an amplification product obtained by amplification of a calibration sequence with a pair of primers designed to produce a bioagent identifying amplicon.

**[0062]** The term "calibration sequence" refers to a polynucleotide sequence to which a given pair of primers hybridizes for the purpose of producing an internal (i.e: included in the reaction) calibration standard amplification product for use in determining the quantity of a bioagent in a sample. The calibration sequence may be expressly added to an amplification reaction, or may already be present in the sample prior to analysis.

**[0063]** The term "clade primer pair" refers to a primer pair designed to produce bioagent identifying amplicons for

species belonging to a clade group. A clade primer pair may also be considered as a speciating primer pair.

**[0064]** The term "codon" refers to a set of three adjoined nucleotides (triplet) that codes for an amino acid or a termination signal.

**[0065]** In context of this invention, the term "codon base composition analysis," refers to determination of the base composition of an individual codon by obtaining a bioagent identifying amplicon that includes the codon. The bioagent identifying amplicon will at least include regions of the target nucleic acid sequence to which the primers hybridize for generation of the bioagent identifying amplicon as well as the codon being analyzed, located between the two primer hybridization regions.

[0066] As used herein, the terms "complementary" or "complementarity" are used in reference to polynucleotides (i.e., a sequence of nucleotides such as an oligonucleotide or a target nucleic acid) related by the base-pairing rules. For example, for the sequence "5'-A-G-T-3'," is complementary to the sequence "3'-T-C-A-5'." Complementarity may be "partial," in which only some of the nucleic acids' bases are matched according to the base pairing rules. Or, there may be "complete" or "total" complementarity between the nucleic acids. The degree of complementarity between nucleic acid strands has significant effects on the efficiency and strength of hybridization between nucleic acid strands. This is of particular importance in amplification reactions, as well as detection methods that depend upon binding between nucleic acids. Either term may also be used in reference to individual nucleotides, especially within the context of polynucleotides. For example, a particular nucleotide within an oligonucleotide may be noted for its complementarity, or lack thereof, to a nucleotide within another nucleic acid strand, in contrast or comparison to the complementarity between the rest of the oligonucleotide and the nucleic acid strand.

[0067] The term "complement of a nucleic acid sequence" as used herein refers to an oligonucleotide which, when aligned with the nucleic acid sequence such that the 5' end of one sequence is paired with the 3' end of the other, is in "antiparallel association." Certain bases not commonly found in natural nucleic acids may be included in the nucleic acids of the present invention and include, for example, inosine and 7-deazaguanine. Complementarity need not be perfect; stable duplexes may contain mismatched base pairs or unmatched bases. Those skilled in the art of nucleic acid technology can determine duplex stability empirically considering a number of variables including, for example, the length of the oligonucleotide, base composition and sequence of the oligonucleotide, ionic strength and incidence of mismatched base pairs. Where a first oligonucleotide is complementary to a region of a target nucleic acid and a second oligonucleotide has complementary to the same region (or a portion of this region) a "region of overlap" exists along the target nucleic acid. The degree of overlap will vary depending upon the extent of the complementarity

**[0068]** In context of this invention, the term "divisionwide primer pair" refers to a primer pair designed to produce bioagent identifying amplicons within sections of a broad spectrum of bioagents For example, primer pair number 367, a division-wide primer pair, is designed to produce bioagent identifying amplicons for the beta-proteobacteria division of bacteria.

**[0069]** As used herein, the term "concurrently amplifying" used with respect to more than one amplification reaction refers to the act of simultaneously amplifying more than one nucleic acid in a single reaction mixture.

**[0070]** As used herein, the term "drill down primer pair" refers to a primer pair designed to produce bioagent identifying amplicons for identification of sub-species characteristics.

[0071] The term "duplex" refers to the state of nucleic acids in which the base portions of the nucleotides on one strand are bound through hydrogen bonding the their complementary bases arrayed on a second strand. The condition of being in a duplex form reflects on the state of the bases of a nucleic acid. By virtue of base pairing, the strands of nucleic acid also generally assume the tertiary structure of a double helix, having a major and a minor groove. The assumption of the helical form is implicit in the act of becoming duplexed.

**[0072]** As used herein, the term "etiology" refers to the causes or origins, of diseases or abnormal physiological conditions.

**[0073]** The term "gene" refers to a DNA sequence that comprises control and coding sequences necessary for the production of an RNA having a non-coding function (e.g., a ribosomal or transfer RNA), a polypeptide or a precursor. The RNA or polypeptide can be encoded by a full length coding sequence or by any portion of the coding sequence so long as the desired activity or function is retained.

[0074] The terms "homology,""homologous" and "sequence identity" refer to a degree of identity. There may be partial homology or complete homology. A partially homologous sequence is one that is less than 100% identical to another sequence. Determination of sequence identity is described in the following example: a primer 20 nucleobases in length which is otherwise identical to another 20 nucleobase primer but having two non-identical residues has 18 of 20 identical residues (18/20=0.9 or 90% sequence identity). In another example, a primer 15 nucleobases in length having all residues identical to a 15 nucleobase segment of a primer 20 nucleobases in length would have 15/20=0.75 or 75% sequence identity with the 20 nucleobase primer. In context of the present invention, sequence identity is meant to be properly determined when the query sequence and the subject sequence are both described and aligned in the 5' to 3' direction. Sequence alignment algorithms such as BLAST, will return results in two different alignment orientations. In the Plus/Plus orientation, both the query sequence and the subject sequence are aligned in the 5' to 3' direction. On the other hand, in the Plus/Minus orientation, the query sequence is in the 5' to 3' direction while the subject sequence is in the 3' to 5' direction. It should be understood that with respect to the primers of the present invention, sequence identity is properly determined when the alignment is designated Plus/Plus. Sequence identity may also encompass alternate or modified nucleobases that perform in a functionally similar manner to the regular nucleobases adenine, thymine, guanine and cytosine with respect to hybridization and primer extension in amplification reactions. In a non-limiting example, if the 5-propynyl pyrimidines propyne C and/or propyne T replace one or more C or T residues in one primer which is otherwise identical to another primer in sequence and length, the two primers will have 100% sequence identity with each other. In another non-limiting example, Inosine (I) may be used as a replacement for G or T and effectively hybridize to C, A or U (uracil). Thus, if inosine replaces one or more C, A or U residues in one primer which is otherwise identical to another primer in sequence and length, the two primers will have 100% sequence identity with each other. Other such modified or universal bases may exist which would perform in a functionally similar manner for hybridization and amplification reactions and will be understood to fall within this definition of sequence identity.

**[0075]** As used herein, "housekeeping gene" refers to a gene encoding a protein or RNA involved in basic functions required for survival and reproduction of a bioagent. House-keeping genes include, but are not limited to genes encoding RNA or proteins involved in translation, replication, recombination and repair, transcription, nucleotide metabolism, amino acid metabolism, lipid metabolism, energy generation, uptake, secretion and the like.

[0076] As used herein, the term "hybridization" is used in reference to the pairing of complementary nucleic acids. Hybridization and the strength of hybridization (i.e., the strength of the association between the nucleic acids) is influenced by such factors as the degree of complementary between the nucleic acids, stringency of the conditions involved, and the T<sub>m</sub> of the formed hybrid. "Hybridization" methods involve the annealing of one nucleic acid to another, complementary nucleic acid, i.e., a nucleic acid having a complementary nucleotide sequence. The ability of two polymers of nucleic acid containing complementary sequences to find each other and anneal through base pairing interaction is a well-recognized phenomenon. The initial observations of the "hybridization" process by Marmur and Lane, Proc. Natl. Acad. Sci. USA 46:453 (1960) and Doty et al., Proc. Natl. Acad. Sci. USA 46:461 (1960) have been followed by the refinement of this process into an essential tool of modern biology.

**[0077]** The term "in silico" refers to processes taking place via computer calculations. For example, electronic PCR (ePCR) is a process analogous to ordinary PCR except that it is carried out using nucleic acid sequences and primer pair sequences stored on a computer formatted medium.

**[0078]** As used herein, "intelligent primers" are primers that are designed to bind to highly conserved sequence regions of a bioagent identifying amplicon that flank an intervening variable region and, upon amplification, yield amplification products which ideally provide enough variability to distinguish individual bioagents, and which are amenable to molecular mass analysis. By the term "highly conserved," it is meant that the sequence regions exhibit between about 80-100%, or between about 90-100%, or between about 95-100% identity among all, or at least 70%, at least 80%, at least 90%, at least 95%, or at least 99% of species or strains.

**[0079]** The "ligase chain reaction" (LCR; sometimes referred to as "Ligase Amplification Reaction" (LAR) described by Barany, Proc. Natl. Acad. Sci., 88:189 (1991); Barany, PCR Methods and Applic., 1:5 (1991); and Wu and

Wallace, Genomics 4:560 (1989) has developed into a well-recognized alternative method for amplifying nucleic acids. In LCR, four oligonucleotides, two adjacent oligonucleotides which uniquely hybridize to one strand of target DNA, and a complementary set of adjacent oligonucleotides, that hybridize to the opposite strand are mixed and DNA ligase is added to the mixture. Provided that there is complete complementarity at the junction, ligase will covalently link each set of hybridized molecules. Importantly, in LCR, two probes are ligated together only when they base-pair with sequences in the target sample, without gaps or mismatches. Repeated cycles of denaturation, hybridization and ligation amplify a short segment of DNA. LCR has also been used in combination with PCR to achieve enhanced detection of single-base changes. However, because the four oligonucleotides used in this assay can pair to form two short ligatable fragments, there is the potential for the generation of target-independent background signal. The use of LCR for mutant screening is limited to the examination of specific nucleic acid positions.

[0080] The term "locked nucleic acid" or "LNA" refers to a nucleic acid analogue containing one or more 2'-O, 4'-Cmethylene- $\beta$ -D-riboftiranosyl nucleotide monomers in an RNA mimicking sugar conformation. LNA oligonucleotides display unprecedented hybridization affinity toward complementary single-stranded RNA and complementary single- or double-stranded DNA. LNA oligonucleotides induce A-type (RNA-like) duplex conformations.

**[0081]** As used herein, the term "mass-modifying tag" refers to any modification to a given nucleotide which results in an increase in mass relative to the analogous non-mass modified nucleotide. Mass-modifying tags can include heavy isotopes of one or more elements included in the nucleotide such as carbon-13 for example. Other possible modifications include addition of substituents such as iodine or bromine at the 5 position of the nucleobase for example.

**[0082]** The term "mass spectrometry" refers to measurement of the mass of atoms or molecules. The molecules are first converted to ions, which are separated using electric or magnetic fields according to the ratio of their mass to electric charge. The measured masses are used to identity the molecules.

**[0083]** The term "microorganism" as used herein means an organism too small to be observed with the unaided eye and includes, but is not limited to bacteria, virus, protozoans, fungi; and ciliates.

**[0084]** The term "multi-drug resistant" or multiple-drug resistant" refers to a microorganism which is resistant to more than one of the antibiotics or antimicrobial agents used in the treatment of said microorganism.

**[0085]** The term "multiplex PCR" refers to a PCR reaction where more than one primer set is included in the reaction pool allowing 2 or more different DNA targets to be amplified by PCR in a single reaction tube.

**[0086]** The term "non-template tag" refers to a stretch of at least three guanine or cytosine nucleobases of a primer used to produce a bioagent identifying amplicon which are not complementary to the template. A non-template tag is incorporated into a primer for the purpose of increasing the primer-duplex stability of later cycles of amplification by

incorporation of extra G-C pairs which each have one additional hydrogen bond relative to an A-T pair.

**[0087]** The term "nucleic acid sequence" as used herein refers to the linear composition of the nucleic acid residues A, T, C or G or any modifications thereof, within an oligonucleotide, nucleotide or polynucleotide, and fragments or portions thereof, and to DNA or RNA of genomic or synthetic origin which may be single or double stranded, and represent the sense or antisense strand

**[0088]** As used herein, the term "nucleobase" is synonymous with other terms in use in the art including "nucleotide,""deoxynucleotide,""nucleotide residue,""deoxynucleotide residue,""nucleotide triphosphate (NTP)," or deoxynucleotide triphosphate (dNTP).

**[0089]** The term "nucleotide analog" as used herein refers to modified or non-naturally occurring nucleotides such as 5-propynyl pyrimidines (i.e., 5-propynyl-dTTP and 5-propynyl-dTCP), 7-deaza purines (i.e., 7-deaza-dATP and 7-deaza-dGTP). Nucleotide analogs include base analogs and comprise modified forms of deoxyribonucleotides as well as ribonucleotides.

[0090] The term "oligonucleotide" as used herein is defined as a molecule comprising two or more deoxyribonucleotides or ribonucleotides, preferably at least 5 nucleotides, more preferably at least about 13 to 35 nucleotides. The exact size will depend on many factors, which in turn depend on the ultimate function or use of the oligonucleotide. The oligonucleotide may be generated in any manner, including chemical synthesis, DNA replication, reverse transcription, PCR, or a combination thereof. Because mononucleotides are reacted to make oligonucleotides in a manner such that the 5' phosphate of one mononucleotide pentose ring is attached to the 3' oxygen of its neighbor in one direction via a phosphodiester linkage, an end of an oligonucleotide is referred to as the "5'-end" if its 5' phosphate is not linked to the 3' oxygen of a mononucleotide pentose ring and as the "3'-end" if its 3' oxygen is not linked to a 5' phosphate of a subsequent mononucleotide pentose ring. As used herein, a nucleic acid sequence, even if internal to a larger oligonucleotide, also may be said to have 5' and 3' ends. A first region along a nucleic acid strand is said to be upstream of another region if the 3' end of the first region is before the 5' end of the second region when moving along a strand of nucleic acid in a 5' to 3' direction. All oligonucleotide primers disclosed herein are understood to be presented in the 5' to 3' direction when reading left to right. When two different, non-overlapping oligonucleotides anneal to different regions of the same linear complementary nucleic acid sequence, and the 3' end of one oligonucleotide points towards the 5' end of the other, the former may be called the "upstream" oligonucleotide and the latter the "downstream" oligonucleotide. Similarly, when two overlapping oligonucleotides are hybridized to the same linear complementary nucleic acid sequence, with the first oligonucleotide positioned such that its 5' end is upstream of the 5' end of the second oligonucleotide, and the 340 end of the first oligonucleotide is upstream of the 3' end of the second oligonucleotide, the first oligonucleotide may be called the "upstream" oligonucleotide and the second oligonucleotide may be called the "downstream" oligonucleotide.

**[0091]** In the context of this invention, a "pathogen" is a bioagent which causes a disease or disorder.

**[0092]** As used herein, the terms "PCR product,""PCR fragment," and "amplification product" refer to the resultant mixture of compounds after two or more cycles of the PCR steps of denaturation, annealing and extension are complete. These terms encompass the case where there has been amplification of one or more segments of one or more target sequences.

[0093] The term "peptide nucleic acid" ("PNA") as used herein refers to a molecule comprising bases or base analogs such as would be found in natural nucleic acid, but attached to a peptide backbone rather than the sugar-phosphate backbone typical of nucleic acids. The attachment of the bases to the peptide is such as to allow the bases to base pair with complementary bases of nucleic acid in a manner similar to that of an oligonucleotide. These small molecules, also designated anti gene agents, stop transcript elongation by binding to their complementary strand of nucleic acid (Nielsen, et al. Anticancer Drug Des. 8:53 63).

**[0094]** The term "polymerase" refers to an enzyme having the ability to synthesize a complementary strand of nucleic acid from a starting template nucleic acid strand and free dNTPs.

[0095] As used herein, the term "polymerase chain reaction" ("PCR") refers to the method of K. B. Mullis U.S. Pat. Nos. 4,683,195, 4,683,202, and 4,965,188, hereby incorporated by reference, that describe a method for increasing the concentration of a segment of a target sequence in a mixture of genomic DNA without cloning or purification. This process for amplifying the target sequence consists of introducing a large excess of two oligonucleotide primers to the DNA mixture containing the desired target sequence, followed by a precise sequence of thermal cycling in the presence of a DNA polymerase. The two primers are complementary to their respective strands of the double stranded target sequence. To effect amplification, the mixture is denatured and the primers then annealed to their complementary sequences within the target molecule. Following annealing, the primers are extended with a polymerase so as to form a new pair of complementary strands. The steps of denaturation, primer annealing, and polymerase extension can be repeated many times (i.e., denaturation, annealing and extension constitute one "cycle"; there can be numerous "cycles") to obtain a high concentration of an amplified segment of the desired target sequence. The length of the amplified segment of the desired target sequence is determined by the relative positions of the primers with respect to each other, and therefore, this length is a controllable parameter. By virtue of the repeating aspect of the process, the method is referred to as the "polymerase chain reaction" (hereinafter "PCR"). Because the desired amplified segments of the target sequence become the predominant sequences (in terms of concentration) in the mixture, they are said to be "PCR amplified." With PCR, it is possible to amplify a single copy of a specific target sequence in genomic DNA to a level detectable by several different methodologies (e.g., hybridization with a labeled probe; incorporation of biotinylated primers followed by avidinenzyme conjugate detection; incorporation of 32P-labeled deoxynucleotide triphosphates, such as dCTP or DATP, into the amplified segment). In addition to genomic DNA, any oligonucleotide or polynucleotide sequence can be amplified with the appropriate set of primer molecules. In particular, the amplified segments created by the PCR process itself are, themselves, efficient templates for subsequent PCR amplifications.

**[0096]** The term "polymerization means" or "polymerization agent" refers to any agent capable of facilitating the addition of nucleoside triphosphates to an oligonucleotide. Preferred polymerization means comprise DNA and RNA polymerases.

**[0097]** As used herein, the terms "pair of primers," or "primer pair" are synonymous. A primer pair is used for amplification of a nucleic acid sequence. A pair of primers comprises a forward primer and a reverse primer. The forward primer hybridizes to a sense strand of a target gene sequence to be amplified and primes synthesis of an antisense strand (complementary to the sense strand) using the target sequence as a template. A reverse primer hybridizes to the antisense strand of a target gene sequence to be amplified and primes synthesis of a sense strand of a target sequence as a template. A reverse primer hybridizes to the antisense strand of a target gene sequence to be amplified and primes synthesis of a sense strand (complementary to the antisense strand) using the target sequence as a template.

[0098] The primers are designed to bind to highly conserved sequence regions of a bioagent identifying amplicon that flank an intervening variable region and yield amplification products which ideally provide enough variability to distinguish each individual bioagent, and which are amenable to molecular mass analysis. In some embodiments, the highly conserved sequence regions exhibit between about 80-100%, or between about 90-100%, or between about 95-100% identity, or between about 99-100% identity. The molecular mass of a given amplification product provides a means of identifying the bioagent from which it was obtained, due to the variability of the variable region. Thus design of the primers requires selection of a variable region with appropriate variability to resolve the identity of a given bioagent. Bioagent identifying amplicons are ideally specific to the identity of the bioagent.

**[0099]** Properties of the primers may include any number of properties related to structure including, but not limited to: nucleobase length which may be contiguous (linked together) or non-contiguous (for example, two or more contiguous segments which are joined by a linker or loop moiety), modified or universal nucleobases (used for specific purposes such as for example, increasing hybridization affinity, preventing non-templated adenylation and modifying molecular mass) percent complementarity to a given target sequences.

[0100] Properties of the primers also include functional features including, but not limited to, orientation of hybridization (forward or reverse) relative to a nucleic acid template. The coding or sense strand is the strand to which the forward priming primer hybridizes (forward priming orientation) while the reverse priming primer hybridizes to the non-coding or antisense strand (reverse priming orientation). The functional properties of a given primer pair also include the generic template nucleic acid to which the primer pair hybridizes. For example, identification of bioagents can be accomplished at different levels using primers suited to resolution of each individual level of identification. Broad range survey primers are designed with the objective of identifying a bioagent as a member of a particular division (e.g., an order, family, genus or other such grouping of bioagents above the species level of bioagents). In some embodiments, broad range survey intelligent primers are

capable of identification of bioagents at the species or sub-species level. Other primers may have the functionality of producing bioagent identifying amplicons for members of a given taxonomic genus, clade, species, sub-species or genotype (including genetic variants which may include presence of virulence genes or antibiotic resistance genes or mutations). Additional functional properties of primer pairs include the functionality of performing amplification either singly (single primer pair per amplification reaction vessel) or in a multiplex fashion (multiple primer pairs and multiple amplification reactions within a single reaction vessel).

**[0101]** As used herein, the terms "purified" or "substantially purified" refer to molecules, either nucleic or amino acid sequences, that are removed from their natural environment, isolated or separated, and are at least 60% free, preferably 75% free, and most preferably 90% free from other components with which they are naturally associated. An "isolated polynucleotide" or "isolated oligonucleotide" is therefore a substantially purified polynucleotide.

**[0102]** The term "reverse transcriptase" refers to an enzyme having the ability to transcribe DNA from an RNA template. This enzymatic activity is known as reverse transcriptase activity. Reverse transcriptase activity is desirable in order to obtain DNA from RNA viruses which can then be amplified and analyzed by the methods of the present invention

**[0103]** The term "Ribosomal RNA" or "rRNA" refers to the primary ribonucleic acid constituent of ribosomes. Ribosomes are the protein-manufacturing organelles of cells and exist in the cytoplasm. Ribosomal RNAs are transcribed from the DNA genes encoding them.

[0104] The term "sample" in the present specification and claims is used in its broadest sense. On the one hand it is meant to include a specimen or culture (e.g., microbiological cultures). On the other hand, it is meant to include both biological and environmental samples. A sample may include a specimen of synthetic origin. Biological samples may be animal, including human, fluid, solid (e.g., stool) or tissue, as well as liquid and solid food and feed products and ingredients such as dairy items, vegetables, meat and meat by-products, and waste. Biological samples may be obtained from all of the various families of domestic animals, as well as feral or wild animals, including, but not limited to, such animals as ungulates, bear, fish, lagamorphs, rodents, etc. Environmental samples include environmental material such as surface matter, soil, water and industrial samples, as well as samples obtained from food and dairy processing instrunents, apparatus, equipment, utensils, disposable and nondisposable items. These examples are not to be construed as limiting the sample types applicable to the present invention. The term "source of target nucleic acid" refers to any sample that contains nucleic acids (RNA or DNA). Particularly preferred sources of target nucleic acids are biological samples including, but not limited to blood, saliva, cerebral spinal fluid, pleural fluid, milk, lymph, sputum and semen.

**[0105]** As used herein, the term "sample template" refers to nucleic acid originating from a sample that is analyzed for the presence of "target" (defined below). In contrast, "background template" is used in reference to nucleic acid other than sample template that may or may not be present in a sample. Background template is often a contaminant. It may be the result of carryover, or it may be due to the presence

of nucleic acid contaminants sought to be purified away from the sample. For example, nucleic acids from organisms other than those to be detected may be present as background in a test sample.

**[0106]** A "segment" is defined herein as a region of nucleic acid within a target sequence.

[0107] The "self-sustained sequence replication reaction" (3SR) (Guatelli et al., Proc. Natl. Acad. Sci., 87:1874-1878 [1990], with an erratum at Proc. Natl. Acad. Sci., 87:7797 [1990]) is a transcription-based in vitro amplification system (Kwok et al., Proc. Natl. Acad. Sci., 86:1173-1177 [1989]) that can exponentially amplify RNA sequences at a uniform temperature. The amplified RNA can then be utilized for mutation detection (Fahy et al., PCR Meth. Appl., 1:25-33 [1991]). In this method, an oligonucleotide primer is used to add a phage RNA polymerase promoter to the 5' end of the sequence of interest. In a cocktail of enzymes and substrates that includes a second primer, reverse transcriptase, RNase H, RNA polymerase and ribo- and deoxyribonucleoside triphosphates, the target sequence undergoes repeated rounds of transcription, cDNA synthesis and second-strand synthesis to amplify the area of interest. The use of 3SR to detect mutations is kinetically limited to screening small segments of DNA (e.g., 200-300 base pairs).

**[0108]** As used herein, the term ""sequence alignment"" refers to a listing of multiple DNA or amino acid sequences and aligns them to highlight their similarities. The listings can be made using bioinformatics computer programs.

**[0109]** In context of this invention, the term "speciating primer pair" refers to a primer pair designed to produce a bioagent identifying amplicon with the diagnostic capability of identifying species members of a group of genera or a particular genus of bioagents. Primer pair number 2922, for example, is a speciating primer pair used to identify species members of the bacterial genus Acinetobacter. Primer pair number 352 is a speciating primer pair used to identify species members of the bacterial genera *Streptococcus*, *Enterococcus*, *Staphylococcus* and *Bacillus*.

**[0110]** In context of this invention, the term "species confirmation primer pair" refers to a primer pair designed to produce a bioagent identifying amplicon with the diagnostic capability to unambiguously produce a unique base composition to identify a particular species of bioagent.

**[0111]** As used herein, a "sub-species characteristic" is a genetic characteristic that provides the means to distinguish two members of the same bioagent species. For example, one viral strain could be distinguished from another viral strain of the same species by possessing a genetic change (e.g., for example, a nucleotide deletion, addition or substitution) in one of the viral genes, such as the RNA-dependent RNA polymerase.

**[0112]** As used herein, the term "target," refers to a nucleic acid sequence or structure to be detected or characterized. Thus, the "target" is sought to be sorted out from other nucleic acid sequences and contains a sequence that has at least partial complementarity with an oligonucleotide primer. The target nucleic acid may comprise single- or double-stranded DNA or RNA. A "segment" is defined as a region of nucleic acid within the target sequence.

**[0113]** The term "template" refers to a strand of nucleic acid on which a complementary copy is built from nucleo-

side triphosphates through the activity of a templatedependent nucleic acid polymerase. Within a duplex the template strand is, by convention, depicted and described as the "bottom" strand. Similarly, the non-template strand is often depicted and described as the "top" strand.

**[0114]** As used herein, the term " $T_m$ " is used in reference to the "melting temperature." The melting temperature is the temperature at which a population of double-stranded nucleic acid molecules becomes half dissociated into single strands. Several equations for calculating the  $T_m$  of nucleic acids are well known in the art. As indicated by standard references, a simple estimate of the  $T_m$  value may be calculated by the equation:  $T_m$ =81.5+0.41 (% G+C), when a nucleic acid is in aqueous solution at 1 M NaCl (see e.g., Anderson and Young, Quantitative Filter Hybridization, in Nucleic Acid Hybridization (1985). Other references (e.g., Allawi, H. T. & SantaLucia, J., Jr. Thermodynamics and NMR of internal G.T mismatches in DNA. Biochemistry 36, 10581-94 (1997) include more sophisticated computations which take structural and environmental, as well as sequence characteristics into account for the calculation of  $T_m$ .

[0115] The term "triangulation genotyping analysis" refers to a method of genotyping a bioagent by measurement of molecular masses or base compositions of amplification products, corresponding to bioagent identifying amplicons, obtained by amplification of regions of more than one gene. In this sense, the term "triangulation" refers to a method of establishing the accuracy of information by comparing three or more types of independent points of view bearing on the same findings. Triangulation genotyping analysis carried out with a plurality of triangulation genotyping analysis primers yields a plurality of base compositions that then provide a pattern or "barcode" from which a species type can be assigned. The species type may represent a previously known sub-species or strain, or may be a previously unknown strain having a specific and previously unobserved base composition barcode indicating the existence of a previously unknown genotype.

**[0116]** As used herein, the term "triangulation genotyping analysis primer pair" is a primer pair designed to produce bioagent identifying amplicons for determining species types in a triangulation genotyping analysis.

**[0117]** The employment of more than one bioagent identifying amplicon for identification of a bioagent is herein referred to as "triangulation identification." Triangulation identification is pursued by analyzing a plurality of bioagent identifying amplicons selected within multiple core genes. This process is used to reduce false negative and false positive signals, and enable reconstruction of the origin of hybrid or otherwise engineered bioagents. For example, identification of the three part toxin genes typical of *B. anthracis* (Bowen et al., J. Appl. Microbiol., 1999, 87, 270-278) in the absence of the expected signatures from the *B. anthracis* genome would suggest a genetic engineering event.

**[0118]** In the context of this invention, the term "unknown bioagent" may mean either: (i) a bioagent whose existence is known (such as the well known bacterial species Staphylococcus aureus for example) but which is not known to be in a sample to be analyzed, or (ii) a bioagent whose existence is not known (for example, the SARS coronavirus was unknown prior to Apr. 2003). For example, if the method for

identification of coronaviruses disclosed in commonly owned U.S. patent Ser. No. 10/829,826 (incorporated herein by reference in its entirety) was to be employed prior to April 2003 to identify the SARS coronavirus in a clinical sample, both meanings of "unknown" bioagent are applicable since the SARS coronavirus was unknown to science prior to April, 2003 and since it was not known what bioagent (in this case a coronavirus) was present in the sample. On the other hand, if the method of U.S. patent Ser. No. 10/829,826 was to be employed subsequent to April 2003 to identify the SARS coronavirus in a clinical sample, only the first meaning (i) of "unknown" bioagent would apply since the SARS coronavirus became known to science subsequent to April 2003 and since it was not known what bioagent was present in the sample.

**[0119]** The term "variable sequence" as used herein refers to differences in nucleic acid sequence between two nucleic acids. For example, the genes of two different bacterial species may vary in sequence by the presence of single base substitutions and/or deletions or insertions of one or more nucleotides. These two forms of the structural gene are said to vary in sequence from one another. In the context of the present invention, "viral nucleic acid" includes, but is not limited to, DNA, RNA, or DNA that has been obtained from viral RNA, such as, for example, by performing a reverse transcription reaction. Viral RNA can either be single-stranded (of positive or negative polarity) or double-stranded.

**[0120]** The term "virus" refers to obligate, ultramicroscopic, parasites incapable of autonomous replication (i.e., replication requires the use of the host cell's machinery). Viruses can survive outside of a host cell but cannon replicate.

**[0121]** The term "wild-type" refers to a gene or a gene product that has the characteristics of that gene or gene product when isolated from a naturally occurring source. A wild-type gene is that which is most frequently observed in a population and is thus arbitrarily designated the "normal" or "wild-type" form of the gene. In contrast, the term "modified", "mutant" or "polymorphic" refers to a gene or gene product that displays modifications in sequence and or functional properties (i.e., altered characteristics) when compared to the wild-type gene or gene product. It is noted that naturally-occurring mutants can be isolated; these are identified by the fact that they have altered characteristics when compared to the wild-type gene or gene product.

**[0122]** As used herein, a "wobble base" is a variation in a codon found at the third nucleotide position of a DNA triplet. Variations in conserved regions of sequence are often found at the third nucleotide position due to redundancy in the amino acid code.

### DETAILED DESCRIPTION OF EMBODIMENTS

### A. Bioagent Identifying Amplicons

**[0123]** The present invention provides methods for detection and identification of unknown bioagents using bioagent identifying amplicons. Primers are selected to hybridize to conserved sequence regions of nucleic acids derived from a bioagent, and which bracket variable sequence regions to yield a bioagent identifying amplicon, which can be amplified and which is amenable to molecular mass determina-

tion. The molecular mass then provides a means to uniquely identify the bioagent without a requirement for prior knowledge of the possible identity of the bioagent. The molecular mass or corresponding base composition signature of the amplification product is then matched against a database of molecular masses or base composition signatures. A match is obtained when an experimentally-determined molecular mass or base composition of an analyzed amplification product is compared with known molecular masses or base compositions of known bioagent identifying amplicons and the experimentally determined molecular mass or base composition is the same as the molecular mass or base composition of one of the known bioagent identifying amplicons. Alternatively, the experimentally-determined molecular mass or base composition may be within experimental error of the molecular mass or base composition of a known bioagent identifying amplicon and still be classified as a match. In some cases, the match may also be classified using a probability of match model such as the models described in U.S. Ser. No. 11/073,362, which is commonly owned and incorporated herein by reference in entirety. Furthermore, the method can be applied to rapid parallel multiplex analyses, the results of which can be employed in a triangulation identification strategy. The present method provides rapid throughput and does not require nucleic acid sequencing of the amplified target sequence for bioagent detection and identification.

**[0124]** Despite enormous biological diversity, all forms of life on earth share sets of essential, common features in their genomes. Since genetic data provide the underlying basis for identification of bioagents by the methods of the present invention, it is necessary to select segments of nucleic acids which ideally provide enough variability to distinguish each individual bioagent and whose molecular mass is amenable to molecular mass determination.

**[0125]** Unlike bacterial genomes, which exhibit conversation of numerous genes (i.e. housekeeping genes) across all organisms, viruses do not share a gene that is essential and conserved among all virus families. Therefore, viral identification is achieved within smaller groups of related viruses, such as members of a particular virus family or genus. For example, RNA-dependent RNA polymerase is present in all single-stranded RNA viruses and can be used for broad priming as well as resolution within the virus family.

**[0126]** In some embodiments of the present invention, at least one viral nucleic acid segment is amplified in the process of identifying the bioagent. Thus, the nucleic acid segments that can be amplified by the primers disclosed herein and that provide enough variability to distinguish each individual bioagent and whose molecular masses are amenable to molecular mass determination are herein described as bioagent identifying amplicons.

**[0127]** In some embodiments of the present invention, bioagent identifying amplicons comprise from about 45 to about 200 nucleobases (i.e. from about 45 to about 200 linked nucleosides), although both longer and short regions may be used. One of ordinary skill in the art will appreciate that the invention embodies compounds of 45, 46, 47, 48, 49, 50, 51, 52, 53, 54, 55, 56, 57, 58, 59, 60, 61, 62, 63, 64, 65, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95, 96, 97, 98, 99, 100, 101, 102, 103, 104, 105, 106, 107, 108, 109,

110, 111, 112, 113, 114, 115, 116, 117, 118, 119, 120, 121, 122, 123, 124, 125, 126, 127, 128, 129, 130, 131, 132, 133, 134, 135, 136, 137, 138, 139, 140, 141, 142, 143, 144, 145, 146, 147, 148, 149, 150, 151, 152, 153, 154, 155, 156, 157, 158, 159, 160, 161, 162, 163, 164, 165, 166, 167, 168, 169, 170, 171, 172, 173, 174, 175, 176, 177, 178, 179, 180, 181, 182, 183, 184, 185, 186, 187, 188, 189, 190, 191, 192, 193, 194, 195, 196, 197, 198, 199, and 200 nucleobases in length, or any range therewithin.

**[0128]** It is the combination of the portions of the bioagent nucleic acid segment to which the primers hybridize (hybridization sites) and the variable region between the primer hybridization sites that comprises the bioagent identifying amplicon.

[0129] In some embodiments, bioagent identifying amplicons amenable to molecular mass determination which are produced by the primers described herein are either of a length, size or mass compatible with the particular mode of molecular mass determination or compatible with a means of providing a predictable fragmentation pattern in order to obtain predictable fragments of a length compatible with the particular mode of molecular mass determination. Such means of providing a predictable fragmentation pattern of an amplification product include, but are not limited to, cleavage with restriction enzymes or cleavage primers, for example. Thus, in some embodiments, bioagent identifying amplicons are larger than 200 nucleobases and are amenable to molecular mass determination following restriction digestion. Methods of using restriction enzymes and cleavage primers are well known to those with ordinary skill in the art.

**[0130]** In some embodiments, amplification products corresponding to bioagent identifying amplicons are obtained using the polymerase chain reaction (PCR) that is a routine method to those with ordinary skill in the molecular biology arts. Other amplification methods may be used such as ligase chain reaction (LCR), low-stringency single primer PCR, and multiple strand displacement amplification (MDA). These methods are also known to those with ordinary skill.

#### B. Primers and Primer Pairs

[0131] In some embodiments the primers are designed to bind to conserved sequence regions of a bioagent identifying amplicon that flank an intervening variable region and yield amplification products which provide variability sufficient to distinguish each individual bioagent, and which are amenable to molecular mass analysis. In some embodiments, the highly conserved sequence regions exhibit between about 80-100%, or between about 90-100%, or between about 95-100% identity, or between about 99-100% identity. The molecular mass of a given amplification product provides a means of identifying the bioagent from which it was obtained, due to the variability of the variable region. Thus, design of the primers involves selection of a variable region with sufficient variability to resolve the identity of a given bioagent. In some embodiments, bioagent identifying amplicons are specific to the identity of the bioagent.

**[0132]** In some embodiments, identification of bioagents is accomplished at different levels using primers suited to resolution of each individual level of identification. Broad range survey primers are designed with the objective of identifying a bioagent as a member of a particular division (e.g., an order, family, genus or other such grouping of

bioagents above the species level of bioagents). In some embodiments, broad range survey intelligent primers are capable of identification of bioagents at the species or sub-species level.

**[0133]** In some embodiments, drill-down primers are designed with the objective of identifying a bioagent at the sub-species level (including strains, subtypes, variants and isolates) based on sub-species characteristics. Drill-down intelligent primers are not always required for identification at the sub-species level because broad range survey intelligent primers may, in some cases provide sufficient identification objective.

[0134] A representative process flow diagram used for primer selection and validation process is outlined in FIG. 1. For each group of organisms, candidate target sequences are identified (200) from which nucleotide alignments are created (210) and analyzed (220). Primers are then designed by selecting appropriate priming regions (230) to facilitate the selection of candidate primer pairs (240). The primer pairs are then subjected to in silico analysis by electronic PCR (ePCR) (300) wherein bioagent identifying amplicons are obtained from sequence databases such as GenBank or other sequence collections (310) and checked for specificity in silico (320). Bioagent identifying amplicons obtained from GenBank sequences (310) can also be analyzed by a probability model which predicts the capability of a given amplicon to identify unknown bioagents such that the base compositions of amplicons with favorable probability scores are then stored in a base composition database (325). Alternatively, base compositions of the bioagent identifying amplicons obtained from the primers and GenBank sequences can be directly entered into the base composition database (330). Candidate primer pairs (240) are validated by testing their ability to hybridize to target nucleic acid by an in vitro amplification by a method such as PCR analysis (400) of nucleic acid from a collection of organisms (410). Amplification products thus obtained are analyzed by gel electrophoresis or by mass spectrometry to confirm the sensitivity, specificity and reproducibility of the primers used to obtain the amplification products (420).

**[0135]** Many of the important pathogens, including the organisms of greatest concern as biowarfare agents, have been completely sequenced. This effort has greatly facilitated the design of primers for the detection of unknown bioagents. The combination of broad-range priming with division-wide and drill-down priming has been used very successfully in several applications of the technology, including environmental surveillance for biowarfare threat agents and clinical sample analysis for medically important pathogens.

**[0136]** Synthesis of primers is well known and routine in the art. The primers may be conveniently and routinely made through the well-known technique of solid phase synthesis. Equipment for such synthesis is sold by several vendors including, for example, Applied Biosystems (Foster City, Calif.). Any other means for such synthesis known in the art may additionally or alternatively be employed.

**[0137]** In some embodiments primers are employed as compositions for use in methods for identification of viral bioagents as follows: a primer pair composition is contacted with nucleic acid (such as, for example, DNA from a DNA

virus, or DNA reverse transcribed from the RNA of an RNA virus) of an unknown viral bioagent. The nucleic acid is then amplified by a nucleic acid amplification technique, such as PCR for example, to obtain an amplification product that represents a bioagent identifying amplicon. The molecular mass of each strand of the double-stranded amplification product is determined by a molecular mass measurement technique such as mass spectrometry for example, wherein the two strands of the double-stranded amplification product are separated during the ionization process. In some embodiments, the mass spectrometry is electrospray Fourier transform ion cyclotron resonance mass spectrometry (ESI-FTICR-MS) or electrospray time of flight mass spectrometry (ESI-TOF-MS). A list of possible base compositions can be generated for the molecular mass value obtained for each strand and the choice of the correct base composition from the list is facilitated by matching the base composition of one strand with a complementary base composition of the other strand. The molecular mass or base composition thus determined is then compared with a database of molecular masses or base compositions of analogous bioagent identifying amplicons for known viral bioagents. A match between the molecular mass or base composition of the amplification product and the molecular mass or base composition of an analogous bioagent identifying amplicon for a known viral bioagent indicates the identity of the unknown bioagent. In some embodiments, the primer pair used is one of the primer pairs of Table 3. In some embodiments, the method is repeated using a different primer pair to resolve possible ambiguities in the identification process or to improve the confidence level for the identification assignment.

**[0138]** In some embodiments, a bioagent identifying amplicon may be produced using only a single primer (either the forward or reverse primer of any given primer pair), provided an appropriate amplification method is chosen, such as, for example, low stringency single primer PCR (LSSP-PCR). Adaptation of this amplification method in order to produce bioagent identifying amplicons can be accomplished by one with ordinary skill in the art without undue experimentation.

**[0139]** In some embodiments, the oligonucleotide primers are broad range survey primers which hybridize to conserved regions of nucleic acid encoding the PB 1 gene or the NUC gene, gene of all (or between 80% and 100%, between 85% and 100%, between 90% and 100% or between 95% and 100%) known adventitious contaminant viruses and produce bioagent identifying amplicons.

**[0140]** In some cases, the molecular mass or base composition of a viral bioagent identifying amplicon defined by a broad range survey primer pair does not provide enough resolution to unambiguously identify a viral bioagent at or below the species level. These cases benefit from further analysis of one or more viral bioagent identifying amplicons generated from at least one additional broad range survey primer pair or from at least one additional division-wide primer pair. The employment of more than one bioagent identifying amplicon for identification of a bioagent is herein referred to as triangulation identification.

**[0141]** In other embodiments, the oligonucleotide primers are division-wide primers which hybridize to nucleic acid encoding genes of species within a genus of viruses. In other

embodiments, the oligonucleotide primers are drill-down primers which enable the identification of sub-species characteristics. Drill down primers provide the functionality of producing bioagent identifying amplicons for drill-down analyses such as strain typing when contacted with nucleic acid under amplification conditions. Identification of such sub-species characteristics is often critical for determining proper clinical treatment of viral infections. In some embodiments, sub-species characteristics are identified using only broad range survey primers and division-wide and drill-down primers are not used.

**[0142]** In some embodiments, the primers used for amplification hybridize to and amplify genomic DNA, DNA of bacterial plasmids, DNA of DNA viruses or DNA reverse transcribed from RNA of an RNA virus.

[0143] In some embodiments, the primers used for amplification hybridize directly to viral RNA and act as reverse transcription primers for obtaining DNA from direct amplification of viral RNA. Methods of amplifying RNA to produce cDNA using reverse transcriptase are well known to those with ordinary skill in the art and can be routinely established without undue experimentation. [1431 In some embodiments, various computer software programs may be used to aid in design of primers for amplification reactions such as Primer Premier 5 (Premier Biosoft, Palo Alto, Calif.) or OLIGO Primer Analysis Software (Molecular Biology Insights, Cascade, Colo.). These programs allow the user to input desired hybridization conditions such as melting temperature of a primer-template duplex for example. In some embodiments, an in silico PCR search algorithm, such as (ePCR) is used to analyze primer specificity across a plurality of template sequences which can be readily obtained from public sequence databases such as GenBank for example. An existing RNA structure search algorithm (Macke et al., Nucl. Acids Res., 2001, 29, 4724-4735, which is incorporated herein by reference in its entirety) has been modified to include PCR parameters such as hybridization conditions, mismatches, and thermodynamic calculations (SantaLucia, Proc. Natl. Acad. Sci. U.S.A., 1998, 95, 1460-1465, which is incorporated herein by reference in its entirety). This also provides information on primer specificity of the selected primer pairs. In some embodiments, the hybridization conditions applied to the algorithm can limit the results of primer specificity obtained from the algorithm. In some embodiments, the melting temperature threshold for the primer template duplex is specified to be 35° C. or a higher temperature. In some embodiments the number of acceptable mismatches is specified to be seven mismatches or less. In some embodiments, the buffer components and concentrations and primer concentrations may be specified and incorporated into the algorithm, for example, an appropriate primer concentration is about 250 nM and appropriate buffer components are 50 mM sodium or potassium and 1.5 mM Mg<sup>2+</sup>.

**[0144]** One with ordinary skill in the art of design of amplification primers will recognize that a given primer need not hybridize with 100% complementarity in order to effectively prime the synthesis of a complementary nucleic acid strand in an amplification reaction. Moreover, a primer may hybridize over one or more segments such that intervening or adjacent segments are not involved in the hybridization event. (e.g., for example, a loop structure or a hairpin structure). The primers of the present invention may com-

prise at least 70%, at least 75%, at least 80%, at least 85%, at least 90%, at least 95% or at least 99% sequence identity with any of the primers listed in Table 3. Thus, in some embodiments of the present invention, an extent of variation of 70% to 100%, or any range therewithin, of the sequence identity is possible relative to the specific primer sequences disclosed herein. Determination of sequence identity is described in the following example: a primer 20 nucleobases in length which is identical to another 20 nucleobase primer having two non-identical residues has 18 of 20 identical residues ( $1^{18}_{20}$ =0.9 or 90% sequence identity). In another example, a primer 15 nucleobase segment of primer 20 nucleobases in length would have  $1^{5}_{20}$ =0.75 or 75% sequence identity with the 20 nucleobase primer.

**[0145]** Percent homology, sequence identity or complementarity, can be determined by, for example, the Gap program (Wisconsin Sequence Analysis Package, Version 8 for UNIX, Genetics Computer Group, University Research Park, Madison Wis.), using default settings, which uses the algorithm of Smith and Waterman (Adv. Appl. Math., 1981, 2, 482-489). In some embodiments, complementarity of primers with respect to the conserved priming regions of viral nucleic acid is between about 70% and about 75% 80%. In other embodiments, homology, sequence identity or complementarity, is between about 75% and about 80%. In yet other embodiments, homology, sequence identity or complementarity, is at least 85%, at least 90%, at least 92%, at least 94%, at least 95%, at least 96%, at least 97%, at least 98%, at least 99% or is 100%.

**[0146]** In some embodiments, the primers described herein comprise at least 70%, at least 75%, at least 80%, at least 85%, at least 90%, at least 92%, at least 94%, at least 95%, at least 96%, at least 98%, or at least 99%, or 100% (or any range therewithin) sequence identity with the primer sequences specifically disclosed herein.

**[0147]** One with ordinary skill is able to calculate percent sequence identity or percent sequence homology and able to determine, without undue experimentation, the effects of variation of primer sequence identity on the fimction of the primer in its role in priming synthesis of a complementary strand of nucleic acid for production of an amplification product of a corresponding bioagent identifying amplicon.

**[0148]** In one embodiment, the primers are at least 13 nucleobases in length. In another embodiment, the primers are less than 36 nucleobases in length.

[0149] In some embodiments of the present invention, the oligonucleotide primers are 13 to 35 nucleobases in length (13 to 35 linked nucleotide residues). These embodiments comprise oligonucleotide primers 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34 or 35 nucleobases in length, or any range therewithin. The present invention contemplates using both longer and shorter primers. Furthermore, the primers may also be linked to one or more other desired moieties, including, but not limited to, affinity groups, ligands, regions of nucleic acid that are not complementary to the nucleic acid to be amplified, labels, etc. Primers may also form hairpin structures. For example, hairpin primers may be used to amplify short target nucleic acid molecules. The presence of the hairpin may stabilize the amplification complex (see e.g., TAQMAN MicroRNA Assays, Applied Biosystems, Foster City, Calif.). **[0150]** In some embodiments, any oligonucleotide primer pair may have one or both primers with less then 70% sequence homology with a corresponding member of any of the primer pairs of Table 3 if the primer pair has the capability of producing an amplification product corresponding to a bioagent identifying amplicon. In other embodiments, any oligonucleotide primer pair may have one or both primers with a length greater than 35 nucleobases if the primer pair has the capability of producing to a bioagent identifying amplification product corresponding to a bioagent identifying an amplification product corresponding to a bioagent identifying amplicon.

**[0151]** In some embodiments, the function of a given primer may be substituted by a combination of two or more primers segments that hybridize adjacent to each other or that are linked by a nucleic acid loop structure or linker which allows a polymerase to extend the two or more primers in an amplification reaction.

[0152] In some embodiments, the primer pairs used for obtaining bioagent identifying amplicons are the primer pairs of Table 3. In other embodiments, other combinations of primer pairs are possible by combining certain members of the forward primers with certain members of the reverse primers. An example can be seen in Table 3 for three primer pair combinations of forward primer POL\_NC003461\_ 2253\_2279\_F (SEQ ID NO: 45), with the reverse primers POL\_NC003461\_2302\_2329\_R (SEQ ID NO: 315), POL\_NC003461\_2320\_2349\_R, or (SEQ ID NO: 200), POL\_NC003461\_2320\_2352\_R (SEQ ID NO: 326). Arriving at a favorable alternate combination of primers in a primer pair depends upon the properties of the primer pair, most notably the size of the bioagent identifying amplicon that would be produced by the primer pair, which should be between about 45 to about 150 nucleobases in length. Alternatively, a bioagent identifying amplicon longer than 150 nucleobases in length could be cleaved into smaller segments by cleavage reagents such as chemical reagents, or restriction enzymes, for example.

**[0153]** In some embodiments, the primers are configured to amplify nucleic acid of a bioagent to produce amplification products that can be measured by mass spectrometry and from whose molecular masses candidate base compositions can be readily calculated.

**[0154]** In some embodiments, any given primer comprises a modification comprising the addition of a non-templated T residue to the 5' end of the primer (i.e., the added T residue does not necessarily hybridize to the nucleic acid being amplified). The addition of a non-templated T residue has an effect of minimizing the addition of non-templated adenosine residues as a result of the non-specific enzyme activity of Taq polymerase (Magnuson et al., Biotechniques, 1996, 21, 700-709), an occurrence which may lead to ambiguous results arising from molecular mass analysis.

**[0155]** In some embodiments of the present invention, primers may contain one or more universal bases. Because any variation (due to codon wobble in the  $3^{rd}$  position) in the conserved regions among species is likely to occur in the third position of a DNA (or RNA) triplet, oligonucleotide primers can be designed such that the nucleotide corresponding to this position is a base which can bind to more than one nucleotide, referred to herein as a "universal nucleobase." For example, under this "wobble" pairing, inosine (I) binds to U, C or A; guanine (G) binds to U or C, and uridine (U) binds to U or C. Other examples of universal nucleobases include nitroindoles such as 5-nitroindole or 3-nitropyrrole (Loakes et al., Nucleosides and Nucleotides,

1995, 14, 1001-1003), the degenerate nucleotides dP or dK (Hill et al.), an acyclic nucleoside analog containing 5-nitroindazole (Van Aerschot et al., Nucleosides and Nucleotides, 1995, 14, 1053-1056) or the purine analog 1-(2deoxy- $\beta$ -D-ribofuranosyl)-imidazole-4-carboxamide (Sala et al., Nucl. Acids Res., 1996, 24, 3302-3306).

[0156] In some embodiments, to compensate for the somewhat weaker binding by the wobble base, the oligonucleotide primers are designed such that the first and second positions of each triplet are occupied by nucleotide analogs that bind with greater affinity than the unmodified nucleotide. Examples of these analogs include, but are not limited to, 2,6-diaminopurine which binds to thymine, 5-propynyluracil which binds to adenine and 5-propynylcytosine and phenoxazines, including G-clamp, which binds to G. Propynylated pyrimidines are described in U.S. Pat. Nos. 5,645, 985, 5,830,653 and 5,484,908, each of which is commonly owned and incorporated herein by reference in its entirety. Propynylated primers are described in U.S Pre-Grant Publication No. 2003-0170682, which is also commonly owned and incorporated herein by reference in its entirety. Phenoxazines are described in U.S. Pat. Nos. 5,502,177, 5,763,588, and 6,005,096, each of which is incorporated herein by reference in its entirety. G-clamps are described in U.S. Pat. Nos. 6,007,992 and 6,028,183, each of which is incorporated herein by reference in its entirety.

**[0157]** In some embodiments, for broad priming of rapidly evolving RNA viruses, primer hybridization is enhanced using primers containing 5-propynyl deoxy-cytidine and deoxy-thymidine nucleotides. These modified primers offer increased affinity and base pairing selectivity.

**[0158]** In some embodiments, non-template primer tags are used to increase the melting temperature  $(T_m)$  of a primer-template duplex in order to improve amplification efficiency. A non-template tag is at least three consecutive A or T nucleotide residues on a primer which are not complementary to the template. In any given non-template tag, A can be replaced by C or G and T can also be replaced by C or G. Although Watson-Crick hybridization is not expected to occur for a non-template tag relative to the template, the extra hydrogen bond in a G-C pair relative to an A-T pair confers increased stability of the primer-template duplex and improves amplification efficiency for subsequent cycles of amplification when the primers hybridize to strands synthesized in previous cycles.

**[0159]** In other embodiments, propynylated tags may be used in a manner similar to that of the non-template tag, wherein two or more 5-propynylcytidine or 5-propynyluridine residues replace template matching residues on a primer. In other embodiments, a primer contains a modified internucleoside linkage such as a phosphorothioate linkage, for example.

**[0160]** In some embodiments, the primers contain massmodifying tags. Reducing the total number of possible base compositions of a nucleic acid of specific molecular weight provides a means of avoiding a persistent source of ambiguity in determination of base composition of amplification products. Addition of mass-modifying tags to certain nucleobases of a given primer will result in simplification of de novo determination of base composition of a given bioagent identifying amplicon from its molecular mass.

**[0161]** In some embodiments of the present invention, the mass modified nucleobase comprises one or more of the following: for example, 7-deaza-2'-deoxyadenosine-5-triph-

osphate, 5-iodo-2'-deoxyuridine-5'-triphosphate, 5-bromo-2'-deoxyuridine-5'-triphosphate, 5-bromo-2'-deoxycytidine-5'-triphosphate, 5-iodo-2'-deoxycytidine-5'-triphosphate, 5-hydroxy-2'-deoxyuridine-5'-triphosphate, 4-thiothymidine-5'-triphosphate, 5-aza-2'-deoxyuridine-5'-triphosphate, 5-fluoro-2'-deoxyuridine-5'-triphosphate, O6-methyl-2'deoxyguanosine-5'-triphosphate, N2-methyl-2'-deoxyguanosine-5'-triphosphate, 8-oxo-2'-deoxyguanosine-5'-triphosphate. In some embodiments, the mass-modified nucleobase comprises <sup>15</sup>N or <sup>13</sup>C or both <sup>15</sup>N and <sup>13</sup>C.

[0162] In some embodiments, multiplex amplification is performed where multiple bioagent identifying amplicons are amplified with a plurality of primer pairs. The advantages of multiplexing are that fewer reaction containers (for example, wells of a 96- or 384-well plate) are needed for each molecular mass measurement, providing time, resource and cost savings because additional bioagent identification data can be obtained within a single analysis. Multiplex amplification methods are well known to those with ordinary skill and can be developed without undue experimentation. However, in some embodiments, one useful and non-obvious step in selecting a plurality candidate bioagent identifying amplicons for multiplex amplification is to ensure that each strand of each amplification product will be sufficiently different in molecular mass that mass spectral signals will not overlap and lead to ambiguous analysis results. In some embodiments, a 10 Da difference in mass of two strands of one or more amplification products is sufficient to avoid overlap of mass spectral peaks.

**[0163]** In some embodiments, as an alternative to multiplex amplification, single amplification reactions can be pooled before analysis by mass spectrometry. In these embodiments, as for multiplex amplification embodiments, it is useful to select a plurality of candidate bioagent identifying amplicons to ensure that each strand of each amplification product will be sufficiently different in molecular mass that mass spectral signals will not overlap and lead to ambiguous analysis results.

C Determination of Molecular Mass of Bioagent Identifying Amplicons

**[0164]** In some embodiments, the molecular mass of a given bioagent identifying amplicon is determined by mass spectrometry. Mass spectrometry has several advantages, not the least of which is high bandwidth characterized by the ability to separate (and isolate) many molecular peaks across a broad range of mass to charge ratio (m/z). Thus mass spectrometry is intrinsically a parallel detection scheme without the need for radioactive or fluorescent labels, since every amplification product is identified by its molecular mass. The current state of the art in mass spectrometry is such that less than femtomole quantities of material can be readily analyzed to afford information about the molecular

contents of the sample. An accurate assessment of the molecular mass of the material can be quickly obtained, irrespective of whether the molecular weight of the sample is several hundred, or in excess of one hundred thousand atomic mass units (amu) or Daltons.

[0165] In some embodiments, intact molecular ions are generated from amplification products using one of a variety of ionization techniques to convert the sample to gas phase. These ionization methods include, but are not limited to, electrospray ionization (ES), matrix-assisted laser desorption ionization (MALDI) and fast atom bombardment (FAB). Upon ionization, several peaks are observed from one sample due to the formation of ions with different charges. Averaging the multiple readings of molecular mass obtained from a single mass spectrum affords an estimate of molecular mass of the bioagent identifying amplicon. Electrospray ionization mass spectrometry (ESI-MS) is particularly useful for very high molecular weight polymers such as proteins and nucleic acids having molecular weights greater than 10 kDa, since it yields a distribution of multiplycharged molecules of the sample without causing a significant amount of fragmentation.

**[0166]** The mass detectors used in the methods of the present invention include, but are not limited to, Fourier transform ion cyclotron resonance mass spectrometry (FT-ICR-MS), time of flight (TOF), ion trap, quadrupole, magnetic sector, Q-TOF, and triple quadrupole.

D. Base Compositions of Bioagent Identifying Amplicons

[0167] Although the molecular mass of amplification products obtained using intelligent primers provides a means for identification of bioagents, conversion of molecular mass data to a base composition signature is useful for certain analyses. As used herein, "base composition" is the exact number of each nucleobase (A, T, C and G) determined from the molecular mass of a bioagent identifying amplicon. In some embodiments, a base composition provides an index of a specific organism. Base compositions can be calculated from known sequences of known bioagent identifying amplicons and can be experimentally determined by measuring the molecular mass of a given bioagent identifying amplicon, followed by determination of all possible base compositions which are consistent with the measured molecular mass within acceptable experimental error. The following example illustrates determination of base composition from an experimentally obtained molecular mass of a 46-mer amplification product originating at position 1337 of the 16S rRNA of Bacillus anthracis. The forward and reverse strands of the amplification product have measured molecular masses of 14208 and 14079 Da, respectively. The possible base compositions derived from the molecular masses of the forward and reverse strands for the B. anthracis products are listed in Table 1.

TABLE 1

-	Possible Base	Compositions for B	. anthracis 46mer	Amplification	n Product
Calc. Mass	Mass Error	Base	Calc. Mass	Mass Error	Base
Forward	Forward	Composition of	Reverse	Reverse	Composition of
Strand	Strand	Forward Strand	Strand	Strand	Reverse Strand
14208.2935	0.056980	A1 G17 C10 T18	14079.2624	0.080600	A0 G14 C13 T19
14208.3160		A1 G20 C15 T10	14079.2849	0.058060	A0 G17 C18 T11
14208.3386		A1 G23 C20 T2	14079.3075	0.035520	A0 G20 C23 T3

TABLE 1-continued

<u> </u>	Possible Base	Compositions for B.	anthracis 46mer	Amplification	n Product
Calc. Mass Forward Strand	Mass Error Forward Strand	Base Composition of Forward Strand	Calc. Mass Reverse Strand	Mass Error Reverse Strand	Base Composition of Reverse Strand
14208.3074 14208.3300 14208.3525 14208.3751 14208.3439 14208.3655 14208.3890 14208.4116 14208.4030 14208.4255 14208.4255 14208.4395 14208.4395	0.065560 0.043020 0.020480 0.029060 0.029060 0.016020 0.038560 0.029980 0.052520 0.075060 0.066480 0.066480 0.066480	A6 G11 C3 T26 A6 G14 C8 T18 A6 G17 C13 T10 A6 G20 C18 T2 A11 G8 C1 T26 A11 G11 C6 T18 A11 G17 C16 T2 A16 G8 C4 T18 A16 G11 C9 T10 A16 G14 C14 T2 A21 G5 C2 T18 A21 G8 C7 T10 —	14079.2538 14079.2764 14079.2989 14079.3214 14079.3129 14079.3354 <b>14079.3579</b> 14079.3355 <b>14079.3579</b> 14079.3805 14079.3494 14079.3719 14079.3494	$\begin{array}{c} 0.089180\\ 0.066640\\ 0.044100\\ 0.021560\\ 0.000980\\ 0.030140\\ 0.007600\\ 0.014940\\ 0.037480\\ 0.006360\\ 0.028900\\ 0.051440\\ 0.073980\\ 0.065400 \end{array}$	A5 G5 C1 T35 A5 G8 C6 T27 A5 G11 C11 T19 A5 G14 C16 T11 A5 G17 C21 T3 A10 G5 C4 T27 A10 G8 C9 T19 A10 G11 C14 T11 A10 G14 C19 T3 A15 G2 C2 T27 A15 G5 C7 T19 A15 G8 C12 T11 A15 G11 C17 T3 A20 G2 C5 T19
	_	_	14079.4309	0.087940	A20 G5 C10 T13

**[0168]** Among the 16 possible base compositions for the forward strand and the 18 possible base compositions for the reverse strand that were calculated, only one pair (shown in bold) are complementary base compositions, which indicates the true base composition of the amplification product. It should be recognized that this logic is applicable for determination of base compositions of any bioagent identifying amplicon, regardless of the class of bioagent from which the corresponding amplification product was obtained.

**[0169]** In some embodiments, assignment of previously unobserved base compositions (also known as "true unknown base compositions") to a given phylogeny can be accomplished via the use of pattern classifier model algorithms. Base compositions, like sequences, vary slightly from strain to strain within species, for example. In some embodiments, the pattern classifier model is the mutational probability model. On other embodiments, the pattern classifier is the polytope model. The mutational probability model and polytope model are both commonly owned and described in U.S. patent application Ser. No. 11/073,362 which is incorporated herein by reference in entirety.

**[0170]** In one embodiment, it is possible to manage this diversity by building "base composition probability clouds" around the composition constraints for each species. This permits identification of organisms in a fashion similar to sequence analysis. A "pseudo four-dimensional plot" can be used to visualize the concept of base composition probability clouds. Optimal primer design requires optimal choice of bioagent identifying amplicons and maximizes the separation between the base composition signatures of individual bioagents. Areas where clouds overlap indicate regions that may result in a misclassification, a problem which is overcome by a triangulation identification process using bioagent identifying amplicons not affected by overlap of base composition probability clouds.

**[0171]** In some embodiments, base composition probability clouds provide the means for screening potential primer pairs in order to avoid potential misclassifications of base compositions. In other embodiments, base composition probability clouds provide the means for predicting the identity of a bioagent whose assigned base composition was not previously observed and/or indexed in a bioagent identifying amplicon base composition database due to evolutionary transitions in its nucleic acid sequence. Thus, in contrast to probe-based techniques, mass spectrometry determination of base composition does not require prior knowledge of the composition or sequence in order to make the measurement.

**[0172]** The present invention provides bioagent classifying information similar to DNA sequencing and phylogenetic analysis at a level sufficient to identify a given bioagent. Furthermore, the process of determination of a previously unknown base composition for a given bioagent (for example, in a case where sequence information is unavailable) has downstream utility by providing additional bioagent indexing information with which to populate base composition databases. The process of future bioagent identification is thus greatly improved as more BCS indexes become available in base composition databases.

#### E. Triangulation Identification

[0173] In some cases, a molecular mass of a single bioagent identifying amplicon alone does not provide enough resolution to unambiguously identify a given bioagent. The employment of more than one bioagent identifying amplicon for identification of a bioagent is herein referred to as "triangulation identification." Triangulation identification is pursued by determining the molecular masses of a plurality of bioagent identifying amplicons selected within a plurality of housekeeping genes. This process is used to reduce false negative and false positive signals, and enable reconstruction of the origin of hybrid or otherwise engineered bioagents. For example, identification of the three part toxin genes typical of B. anthracis (Bowen et al., J. Appl. Microbiol., 1999, 87, 270-278) in the absence of the expected signatures from the B. anthracis genome would suggest a genetic engineering event.

**[0174]** In some embodiments, the triangulation identification process can be pursued by characterization of bioagent identifying amplicons in a massively parallel fashion using the polymerase chain reaction (PCR), such as multiplex PCR where multiple primers are employed in the same amplification reaction mixture, or PCR in multi-well plate format wherein a different and unique pair of primers is used in multiple wells containing otherwise identical reaction mixtures. Such multiplex and multi-well PCR methods are well known to those with ordinary skill in the arts of rapid throughput amplification of nucleic acids. In other related embodiments, one PCR reaction per well or container may be carried out, followed by an amplicon pooling step wherein the amplification products of different wells are combined in a single well or container which is then subjected to molecular mass analysis. The combination of pooled amplicons can be chosen such that the expected ranges of molecular masses of individual amplicons are not overlapping and thus will not complicate identification of signals.

F. Codon Base Composition Analysis

**[0175]** In some embodiments of the present invention, one or more nucleotide substitutions within a codon of a gene of an infectious organism confer drug resistance upon an organism which can be determined by codon base composition analysis. The organism can be a bacterium, virus, fungus or protozoan.

**[0176]** In some embodiments, the amplification product containing the codon being analyzed is of a length of about 35 to about 150 nucleobases. The primers employed in obtaining the amplification product can hybridize to upstream and downstream sequences directly adjacent to the codon, or can hybridize to upstream and downstream sequence positions away from the codon. The primers may have between about 70% to 100% sequence complementarity with the sequence of the gene containing the codon being analyzed.

**[0177]** In some embodiments, the codon base composition analysis is undertaken

**[0178]** In some embodiments, the codon analysis is undertaken for the purpose of investigating genetic disease in an individual. In other embodiments, the codon analysis is undertaken for the purpose of investigating a drug resistance mutation or any other deleterious mutation in an infectious organism such as a bacterium, virus, fungus or protozoan. In some embodiments, the virus is an adventitious virus identified in a biological product.

**[0179]** In some embodiments, the molecular mass of an amplification product containing the codon being analyzed is measured by mass spectrometry. The mass spectrometry or matrix-assisted laser desorption ionization (MALDI) mass spectrometry. Time-of-flight (TOF) is an example of one mode of mass spectrometry compatible with the analyses of the present invention.

**[0180]** The methods of the present invention can also be employed to determine the relative abundance of drug resistant strains of the organism being analyzed. Relative abundances can be calculated from amplitudes of mass spectral signals with relation to internal calibrants. In some embodiments, known quantities of internal amplification calibrants can be included in the amplification reactions and abundances of analyte amplification product estimated in relation to the known quantities of the calibrants. **[0181]** In some embodiments, upon identification of one or more drug-resistant strains of an infectious organism infecting an individual, one or more alternative treatments can be devised to treat the individual.

G. Determination of the Quantity of a Bioagent

[0182] In some embodiments, the identity and quantity of an unknown bioagent can be determined using the process illustrated in FIG. 2. Primers (500) and a known quantity of a calibration polynucleotide (505) are added to a sample containing nucleic acid of an unknown bioagent. The total nucleic acid in the sample is then subjected to an amplification reaction (510) to obtain amplification products. The molecular masses of amplification products are determined (515) from which are obtained molecular mass and abundance data. The molecular mass of the bioagent identifying amplicon (520) provides the means for its identification (525) and the molecular mass of the calibration amplicon obtained from the calibration polynucleotide (530) provides the means for its identification (535). The abundance data of the bioagent identifying amplicon is recorded (540) and the abundance data for the calibration data is recorded (545), both of which are used in a calculation (550) which determines the quantity of unknown bioagent in the sample.

[0183] A sample comprising an unknown bioagent is contacted with a pair of primers that provide the means for amplification of nucleic acid from the bioagent, and a known quantity of a polynucleotide that comprises a calibration sequence. The nucleic acids of the bioagent and of the calibration sequence are amplified and the rate of amplification is reasonably assumed to be similar for the nucleic acid of the bioagent and of the calibration sequence. The amplification reaction then produces two amplification products: a bioagent identifying amplicon and a calibration amplicon. The bioagent identifying amplicon and the calibration amplicon should be distinguishable by molecular mass while being amplified at essentially the same rate. Effecting differential molecular masses can be accomplished by choosing as a calibration sequence, a representative bioagent identifying amplicon (from a specific species of bioagent) and performing, for example, a 2-8 nucleobase deletion or insertion within the variable region between the two priming sites. The amplified sample containing the bioagent identifying amplicon and the calibration amplicon is then subjected to molecular mass analysis by mass spectrometry, for example. The resulting molecular mass analysis of the nucleic acid of the bioagent and of the calibration sequence provides molecular mass data and abundance data for the nucleic acid of the bioagent and of the calibration sequence. The molecular mass data obtained for the nucleic acid of the bioagent enables identification of the unknown bioagent and the abundance data enables calculation of the quantity of the bioagent, based on the knowledge of the quantity of calibration polynucleotide contacted with the sample.

**[0184]** In some embodiments, construction of a standard curve where the amount of calibration polynucleotide spiked into the sample is varied provides additional resolution and improved confidence for the determination of the quantity of bioagent in the sample. The use of standard curves for analytical determination of molecular quantities is well known to one with ordinary skill and can be performed without undue experimentation.

**[0185]** In some embodiments, multiplex amplification is performed where multiple bioagent identifying amplicons are amplified with multiple primer pairs which also amplify the corresponding standard calibration sequences. In this or other embodiments, the standard calibration sequences are optionally included within a single vector which functions as the calibration polynucleotide. Multiplex amplification methods are well known to those with ordinary skill and can be performed without undue experimentation.

**[0186]** In some embodiments, the calibrant polynucleotide is used as an internal positive control to confirm that amplification conditions and subsequent analysis steps are successful in producing a measurable amplicon. Even in the absence of copies of the genome of a bioagent, the calibration polynucleotide should give rise to a calibration amplicon. Failure to produce a measurable calibration amplicon indicates a failure of amplification or subsequent analysis step such as amplicon purification or molecular mass determination. Reaching a conclusion that such failures have occurred is in itself, a useful event.

**[0187]** In some embodiments, the calibration sequence is comprised of DNA. In some embodiments, the calibration sequence is comprised of RNA.

[0188] In some embodiments, the calibration sequence is inserted into a vector that itself functions as the calibration polynucleotide. In some embodiments, more than one calibration sequence is inserted into the vector that functions as the calibration polynucleotide. Such a calibration polynucleotide is herein termed a "combination calibration polynucleotide." The process of inserting polynucleotides into vectors is routine to those skilled in the art and can be accomplished without undue experimentation. Thus, it should be recognized that the calibration method should not be limited to the embodiments described herein. The calibration method can be applied for determination of the quantity of any bioagent identifying amplicon when an appropriate standard calibrant polynucleotide sequence is designed and used. The process of choosing an appropriate vector for insertion of a calibrant is also a routine operation that can be accomplished by one with ordinary skill without undue experimentation.

H. Identification of Adventitious Viruses

**[0189]** In other embodiments of the present invention, the primer pairs produce bioagent identifying amplicons within stable and highly conserved regions of adventitious contaminant viruses. The advantage to characterization of an amplicon in a highly conserved region is that there is a low probability that the region will evolve past the point of primer recognition, in which case, the amplification step would fail. Such a primer set is thus useful as a broad range survey-type primer. In another embodiment of the present invention, the intelligent primers produce bioagent identifying amplicons in a region which evolves more quickly than the stable region described above. The advantage of characterization bioagent identifying amplicon corresponding to an evolving genomic region is that it is useful for distinguishing emerging strain variants.

**[0190]** The present invention also has significant advantages as a platform for identification of diseases caused by emerging viruses. The present invention eliminates the need for prior knowledge of bioagent sequence to generate hybridization probes. Thus, in another embodiment, the present invention provides a means of determining the etiology of a virus infection when the process of identification of viruses is carried out in a clinical setting and, even when the virus is a new species never observed before. This is possible because the methods are not confounded by naturally occurring evolutionary variations (a major concern for characterization of viruses which evolve rapidly) occurring in the sequence acting as the template for production of the bioagent identifying amplicon. Measurement of molecular mass and determination of base composition is accomplished in an unbiased manner without sequence prejudice.

**[0191]** Another embodiment of the present invention also provides a means of tracking the spread of any species or strain of virus when a plurality of samples obtained from different locations are analyzed by the methods described above in an epidemiological setting. In one embodiment, a plurality of samples from a plurality of different locations is analyzed with primer pairs which produce bioagent identifying amplicons, a subset of which contains a specific virus. The corresponding locations of the members of the virus-containing subset indicate the spread of the specific virus to the corresponding locations.

**[0192]** Members of the Parvoviridae family are small single stranded DNA viruses with genomes of about 4-5 kilobases long. They can be divided into (i) Dependovirus genus that includes the human helper-dependent adeno-associated virus (AAV) serotypes 1 to 8 and the autonomous avian parvoviruses; the adeno associated viruses (AAV 1-8); (ii) Erythrovirus genus that includes the bovine, chipmunk, and autonomous primate parvoviruses, including human viruses B19 and V9; and (iii) Parvovirus genus that include parvoviruses of other animals and rodents (except for chipmunks), carnivores, and pigs, including murine minute virus (MMV). These parvoviruses can infect several cell types and have been described in clinical samples. AAVs in particular, have been implicated in decreased replication, propagation, and growth of other virus.

**[0193]** Exogenous retroviruses are known to cause various malignant and non-malignant diseases in animals over a wide range of species. These viruses infect most known animals and rodents. Examples include, but are not limited to: Deltaretroidvirus (HTLV 1-4, STLV 1-3), Gammaretrovirus (Murine leukemia virus, PERV), Alpharetrovirus: (Avian leucosis virus and Avian endogenous virus) and Human immunodeficiency viruses 1 and 2).

**[0194]** Polyomaviruses are small double-stranded DNA viruses that can infect several species including humans, primates, rodents, rabbits and birds. Because of their tumorigenic and oncogenic potential, it is important to test for these viruses in cell substrates used for vaccine production.

**[0195]** The Papillomaviridae family of viruses contains more that 150 known species representing varying hostspecificity and sequence homology. They have been identified in mammals (humans, simians, bovines, canines, ovines) and in birds. Majority of the human Papillomaviruses (HPVs), including all HPV types traditionally called genital and mucosal HPVs belong to supergroup A. Within supergroup A, there are 11 groups; the most medically important of these are the human Papillomaviruses HPV 16, HPV 18, HPV 31, HPV 45, HPV 11, HPV 6 and HPV 2. Each of these has been reported as "high risk" viruses in the medical literature.

[0196] Other viral families which are potential adventitious contaminants include, but are not limited to: Herpesviridae (Human herpesviruses 1 through 8, Bovine herpesvirus, Canine herpesvirus and Simian cytomegalovirus), Hepadnaviridae (Hepatitis B virus), Hepeviridae (Hepatitis E virus), Deltavirus (Hepatitis delta virus), Adenoviridae (Human adenoviruses A-F and murine adenovirus), Flaviviridae (Bovine viral diarrhea virus, TBE, Yellow fever virus, Dengue viruses 14, WNV and hepatitis C virus), Paramyxoviridae (Pneumonia virus of mice, Sendai virus, and Simian parainfluenza virus 5), Togaviridae (Western equine encephalomyelitis virus), Picornaviridae (Polio (types 1-13), Human hepatitis A, Human coxsackievirus, Human cardiovirus, Human rhinovirus and Bovine rhinovirus), Reoviridae (Mouse rotavirus, reovirus type 3 and Colorado tick fever virus), and Rhabdoviridae (vesicular stomatitis virus).

### I. Kits

**[0197]** The present invention also provides kits for carrying out the methods described herein. In some embodiments, the kit may comprise a sufficient quantity of one or more primer pairs to perform an amplification reaction on a target polynucleotide from a bioagent to form a bioagent identifying amplicon. In some embodiments, the kit may comprise from one to fifty primer pairs, from one to twenty primer pairs, from one to ten primer pairs, or from two to five primer pairs. In some embodiments, the kit may comprise one or more primer pairs recited in Table 3.

[0198] In some embodiments, the kit comprises one or more broad range survey primer(s), division wide primer(s), or drill-down primer(s), or any combination thereof. If a given problem involves identification of a specific bioagent, the solution to the problem may require the selection of a particular combination of primers to provide the solution to the problem. A kit may be designed so as to comprise particular primer pairs for identification of a particular bioagent. A drill-down kit may be used, for example, to distinguish different sub-species types of adventitious contaminant viruses or genetically engineered adventitious contaminant viruses. In some embodiments, the primer pair components of any of these kits may be additionally combined to comprise additional combinations of broad range survey primers and division-wide primers so as to be able to identify the adventitious contaminant virus.

**[0199]** In some embodiments, the kit contains standardized calibration polynucleotides for use as internal amplification calibrants. Internal calibrants are described in commonly owned U.S. patent application Ser. No: 60/545,425 which is incorporated herein by reference in its entirety.

**[0200]** In some embodiments, the kit comprises a sufficient quantity of reverse transcriptase (if an RNA virus is to be identified for example), a DNA polymerase, suitable nucleoside triphosphates (including alternative dNTPs such as inosine or modified dNTPs such as the 5-propynyl pyrimidines or any dNTP containing molecular mass-modifing tags such as those described above), a DNA ligase, and/or reaction buffer, or any combination thereof, for the amplification processes described above. A kit may further include instructions pertinent for the particular embodiment of the kit, such instructions for operation of the method. A kit may also comprise amplification reaction containers such as

microcentrifuge tubes and the like. A kit may also comprise reagents or other materials for isolating bioagent nucleic acid or bioagent identifying amplicons from amplification, including, for example, detergents, solvents, or ion exchange resins which may be linked to magnetic beads. A kit may also comprise a table of measured or calculated molecular masses and/or base compositions of bioagents using the primer pairs of the kit.

**[0201]** In some embodiments, the kit includes a computer program stored on a computer formatted medium (such as a compact disk or portable USB disk drive, for example) comprising instructions which direct a processor to analyze data obtained from the use of the primer pairs of the present invention. The instructions of the software transform data related to amplification products into a molecular mass or base composition which is a useful concrete and tangible result used in identification and/or classification of bioagents. In some embodiments, the kits of the present invention contain all of the reagents sufficient to carry out one or more of the methods described herein.

**[0202]** While the present invention has been described with specificity in accordance with certain of its embodiments, the following examples serve only to illustrate the invention and are not intended to limit the same. In order that the invention disclosed herein may be more efficiently understood, examples are provided below. It should be understood that these examples are for illustrative purposes only and are not to be construed as limiting the invention in any manner.

#### **EXAMPLES**

#### Example 1

### Design and Validation of Primers that Define Bioagent Identifying Amplicons for Adventitious Contaminant Viruses

#### A. General Process of Primer Design

**[0203]** For design of primers that define adventitious contaminant virus identifying amplicons, a series of adventitious contaminant virus genome segment sequences were obtained, aligned and scanned for regions where pairs of PCR primers would amplify products of about 45 to about 150 nucleotides in length and distinguish species and/or individual strains from each other by their molecular masses or base compositions. A typical process shown in **FIG. 1** is employed for this type of analysis.

**[0204]** A database of expected base compositions for each primer region was generated using an in silico PCR search algorithm, such as (ePCR). An existing RNA structure search algorithm (Macke et al., Nucl. Acids Res., 2001, 29, 4724-4735, which is incorporated herein by reference in its entirety) has been modified to include PCR parameters such as hybridization conditions, mismatches, and thermodynamic calculations (SantaLucia, Proc. Natl. Acad. Sci. U.S.A., 1998, 95, 1460-1465, which is incorporated herein by reference in its entirety). This also provides information on primer specificity of the selected primer pairs.

B. Design of Primers for Identification of Parvoviruses

**[0205]** Primer pairs were designed which broadly target the various genera/species of the Parvovirinae family. Par-

voviruses of the most medical concern are: B19, AAV-5 and murine minute virus. Approximately 500 complete Parvovirus genome sequences were obtained from GenBank. These genome sequences (each approximately 5 kilobases long) were aligned and scanned for conserved target regions. Initial survey of the genome alignments revealed very little homology across the three major genera described above. However, regions were identified with significant homologies within each genus that were the target for primer design. In all three genera, the regions of conservation were within two major nodes, one in the rep gene, encoding NS1 protein, and the other in the capsid, cap gene, encoding glycoprotein VP1 protein. This ability to prime across all known instances of species within each of these groups will enable surveillance for known parvoviruses and detection of previously unknown parvoviruses in cell lines.

#### C. Design of Primers for Identification of Retroviruses

[0206] The objective of primer design for this viral family was to obtain primer pairs that would prime and produce retrovirus identifying amplicons for all known members of each of the genus groups and as yet unknown variants. The T-lymphotropic viruses, members of the deltaretrovirus genus, infect primates and cause leukemia and neurologic diseases. These 9 kilobase single stranded RNA viruses are highly transmissible. Primer pairs targeting the transcription activator (tax) gene were designed to broadly prime and resolve all known primate T-lymphotropic viruses including human T-lymphotropic viruses (HTLV-1 and -2 and the newly discovered HTLV-3 and -4), and simian T-lymphotropic viruses (STLV-1, -2 and -3). These primer pairs produce retrovirus identifying amplicons of simian and human T-lymphotropic virus species with distinct base compositions indicating that the primer pairs can yield amplification products which are distinguishable from each other on the basis of molecular masses and base compositions.

#### D. Design of Primers for Identification of Polyomaviruses

[0207] Approximately 200 complete Polyomavirus genome sequences were obtained from GenBank. These genome sequences (approximately 5.3 kilobases long) were aligned to each other using bioinformatics tools built inhouse, and scanned for conserved target regions. Initial survey of the genome alignments revealed a high degree of homology between the primate (SV40) and the human viruses (BK and JC), whereas the rest of the species were highly divergent and did not share much sequence homology with these above species. For primer design purposes, SV40, BK and JC viral species were classified as a group. Nine different primer pairs (primer pairs 2549-2557) were designed to this cluster (See Table 3) and are expected to provide redundant detection and resolution of the three important Polyomavirus species. Most of these primers were targeted to the large T antigen gene of Polyomavirus. Three additional primer pairs (primer pair numbers 2559-2561) were designed to include Lymphotropic papovavirus (LPV, the African green monkey papovavirus). While these new primers were less conserved across any one species, they would nonetheless provide broader coverage of viral detection within this family. Additional primer pairs (RS 10-14) targeting the rest of the viral species (murine, avian, bovine, etc.) were also designed. Taken together, these primers would provide complete coverage of all known Polyomaviruses.

[0208] All of the polyomavirus primer pairs were tested against multiple target species for performance and sensitivity. To test the performance of these primers, plasmid clones containing full length SV40 (ATCC: VRMC-4) and JC virus (ATCC: VRMC-1) DNA were obtained from ATCC. Plasmid concentrations were determined by optical density measurements and used as an approximate estimate of the amount of input viral DNA template. Serial 10-fold dilutions of the plasmid were used for estimating limits of detection. These were tested against the entire panel of 12 primer pairs (primer pair numbers 2549-2561). The primer pairs were initially tested at  $10^{-7}$  and  $10^{-8}$  fold dilutions of each of the plasmids and showed reliable detections, with the exception of primer 2555. Additional testing of a subset of these primer pairs showed that while several primer pairs were able to detect additional, lower dilutions, some of the primer pairs were unsuccessful at producing polyomavirus identifying amplicons below the  $10^{-9}$  dilution. Based on this initial test, a panel of six primers (primer pair numbers 2550, 2551, 2553, 2554, 2557 and 2559) was chosen for use in cell-line characterization. These primers will be tested against known cell lines containing SV-40 and other Polyomaviruses.

**[0209]** For routine screening of cell lines, it is anticipated that as few as two of the primer pairs described above along with the four primers targeting non-human Polyomavirus can provide complete coverage of all known and potentially novel Polyomaviruses. A nucleic acid segment within the large tumor antigen gene provides opportunities for broad priming across human and simian species due to a codon deletion at position 32 of the simian virus 40, which is exemplified by primer pair number 2555 (SEQ ID NO: 112:207). Murine pneumonotropic virus, African green monkey PyV virus, SV40 virus, BK virus, JC virus, hamster PyV and murine PyV virus can be distinguished from each other on the basis of base compositions of amplification products produced with primer pair number 2560 (SEQ ID NOs: 12:260).

#### E. Design of Primers for Identification of Papillomaviruses

**[0210]** Broad primer pairs covering a set of important human Papillomaviruses (HPV 16, 18, 31, 45, 11, 6, 2) were designed (primer pair numbers 2533-2536). These belong to different groups, but have all been reported in literature to be "high risk" Covering all of these species broadly combined with group-specific primer pairs described above would be of great value. Additionally, several primer pairs were designed to cover broadly within a single group or across multiple groups of Papillomaviruses to increase robustness of detection.

**[0211]** All of the primer pairs were tested against a panel of Papillomaviruses obtained from ATCC. The following viruses were obtained as full-length plasmid clones: ATCC 45150D (HPV-6b); ATCC 45151D (HPV-11); ATCC 45152D (HPV-18); and ATCC 45113D (HPV-16). Two of the broad primer pairs (numbers 2534 and 2536) amplified all four viruses tested at two different dilutions of the plasmids. Primer pair number 2535 (SEQ ID NOs: 28:253) amplified only two of the test isolates, while primer pair 2533 (30:268) did not amplify any of the viruses tested.

Based on these initial results, Primer pair numbers 2534 (SEQ ID NOs: 30:267) and 2536 (SEQ ID NOs: 19:267) were selected for further optimization. A series of primer modifications, including, for example, inosine substitutions to overcome potential sequence mismatches were introduced into the forward and reverse primer pairs. Most of the modified primers tested showed improved performance across the test isolates. In addition to the primers broadly targeting the major species, a series of primers targeting Papillomavirus groups, A7, A9 and A10 that account for over 30 different Papillomaviruses were also tested. Table 2 provides the primer pairs used for Papillomavirus identification and indicates isolates tested, target virus groups and major species covered.

TABLE 2

P1	Primer Pairs Targeting Human Papillomaviruses							
Primer Pair Number	Isolates Tested	Target Virus Group	Major Species Covered					
2537 2539 2540	HPV-16	Group A9	HPV-16, HPV-31, HPV-33, HPV-35, HPV-52, HPV-58, HPV-67, and RhPV					
2543 2544 2545	HPV-18	Group A7	HPV-18, HPV-39, HPV-45, HPV-59, HPV-68, and HPV-70					
2546 2547 2548	HPV-6, HPV-11	Group A10	HPV-6, HPV-11, HPV-13, HPV-44, HPV-55, and PCPV					
2541 2542	HPV-6, HPV-11, HPV-18	Groups A1, A7, A8, A10 and A11	>30 different Papilloma viruses					

F. Validation of Primer Pairs Designed for Identification of Papillomaviruses

**[0212]** For additional testing and validation, two different HeLa cell lines infected with HPV-18 were obtained from ATCC (CCL-2 and CCL-2.2). These were tested at limiting dilutions using a subset of the primers listed above. Results are shown below. The primer pairs used for this test included the major human PaV primer pairs, 2534 (SEQ ID NOS: 30:267), 2536 (SEQ ID NOS: 19:267) and 2685 (SEQ ID NOS: 18:272), the multi-group primer 2542 (SEQ ID NOS: 8:294) and 2545 SEQ ID NOS: 98:193) and the Group A10 primer 2546 (SEQ ID NOS: 64:302).

**[0213]** In addition to testing the performance of the primers on the cell lines, plasmid DNA containing HPV-6b was

spiked into the CCL-2 cell line to determine the dynamic range of detection of the two viruses, cell line derived HPV-18 and the plasmid-derived HPV-6b, simultaneously, In all the tests done, the broad primers as well as the Group A7 primers showed detection of HPV-1 8 in both cell lines at input levels between 1-10 cells per well. At an estimated copy number of approximately 20 HPV-18 genomes per cell, this corresponds to detection sensitivities between 20-200 genomes from cell lines containing papillomavirus sequences. In experiments done with a co-spike of HPV-6b plasmid into these cell lines, the detection ranges were comparable. HPV-6b was spiked in at two different, fixed concentrations of 200 copies and 2000 copies per well and amplified with the broad primer pair number 2534. Simultaneous detection of HPV-6b and HPV-18 was observed when the plasmid DNA was spiked in at 2000 copies into a range of CCL-2 cell concentration from 1000 to 0 per well. HPV-18 was detected in all wells with the exception of the lowest input level (10 cells/well), in the presence of 2000 copies of HPV-6b. HPV-6b (2000 copies) was detected in the presence of HeLa cell loads up to 600 cells/well, with an effective HPV-18 concentration of approximately 12000 genomes/well. In another experiment, a plasmid spike of approximately 200 copies per well was used. In this case, HPV-18 was detected at all test concentrations, including the lowest cell concentration of 10 cells per well. The dynamic range for detection of the two viruses simultaneously is between 5-10 fold at the lower and higher ends, giving an overall dynamic range of ~25 fold for the detection of competing templates in the presence of each other. These experiments indicate that two or more viruses can be simultaneously detected using the same assay.

G. Primer Pair Compositions for Identification of Adventitious Viruses

**[0214]** A total of 224 primer pairs were designed. Table 3 which represents a collection of primers (sorted by primer pair number) designed to identify adventitious contaminant viruses using the methods described herein. "I" represents inosine. Tp represents propynylated T and Cp represents propynylated CP, wherein the propynyl substituent is located at the 5-position of the pyrimidine nucleobase. The primer pair number is an in-house database index number. The forward or reverse primer name shown in Table 3 indicates the gene region of the viral genome to which the primer hybridizes relative to a reference sequence. The forward primer name RVL\_X03614\_2256\_2279\_F indicates that the forward primer (\_F) hybridizes to residues 2256-2279 of a respirovirus (Paramyxoviridae) sequence (GenBank Accession r X03614).

TABLE 3

	Pi	rimer Pairs for Identific	cation of Adve	ntitious Contami	nant Viruses	
Primer	~		Forward	3		Reverse
	-			-		
Pair			SEQ ID			SEQ ID
	Forward Primer			Reverse Primer		
Number	Name	Forward Sequence	NO:	Name	Reverse Sequence	NO:
377	RVL_X03614_2256 2279 F	5_ TAGTGCAATCCATCTAGCAGC	CTGT 39	RVL_X03614_2302_ 2324 R	TGATTGTCGCCTTGAACCATT GC	293

Primer Pair	2		Forward SEQ ID			Reverse SEQ ID
Number	Forward Primer Name	Forward Sequence	NO:	Reverse Primer Name	Reverse Sequence	NO:
378	RVL_X03614_2239_ 2261_F	TGGACATTCATCTCTATCAGTGC	132	RVL_X03614_2302_ 2324_R	TGATTGTCGCCTTGAACCATT GC	293
379	PVL_U50363_3176_ 3199_F	TAGATCCACAAGCTTTAGGGTCTG	22	PVL_U50363_3272_ 3296_R	TGTTGTGCACTTTTGGAGAAT ATTT	350
380	PVL_U50363_3153_ 3175_F	TGCTGAATTCGTAACATTGATGA	130	PVL_U50363_3265_ 3286_R	TTTTGGCGAATATTTTGTTTG G	371
381	MVL_AF266286_ 1625_1649_F	TAGGGAGACTTTTTGCTAAAATGAC	33	MVL_AF266286_ 1720_1739_R	TCCTTTGCCATCCCATTGTC	248
382	MVL_AF266286_ 2033_2054_F	TGTTTGCACAGAGGCTAAATGA	170	MVL_AF266286_ 2122_2139_R	TGGGCAATGAGGGTCACT	318
383	MVL_AF266286_ 2002_2024_F	TGCCTTAATTGGAGATATGAGAC	126	MVL_AF266286_ 2035_2054_R	TCATTTAGCCTCTGTGCAAA	234
384	MSVL_AF266286_ 3529_3548_F	TGTGTCATCTGCGAGTGTGG	168	MSVL_AF266286_ 3587_3608_R	TTGTCAATATCATCCAGTTGG C	366
385	PNVL_U50363_ 1285_1301_F	TTTGGACACCCAATGGT	185	PNVL_U50363_ 1318_1335_R	TACTCTAACAGCATCCAT	199
386	PNVL_U50363_ 1878_1898_F	TAGAATGTTTGCTATGCAACC	20	PNVL_U50363_ 1927_1949_R	TTCTCAGCTAACAATTTCTCA GC	357
387	PNVL_U50363_ 1881_1898_F	TATGTTTGCTATGCAACC	51	PNVL_U50363_ 1927_1944_R	TGCTAACAATTTCTCAGC	304
388	PNVL_U50363_ 2669_2687_F	TAGGACCGTGGATAAACAC	27	PNVL_U50363_ 2743_2763_R	TCTTTCCCCTCTGTATTCTAA	278
389	MNVL_NC004148_ 24_47_F	TGTTAATGTCTATCTTCCTGACTC	169	MNVL_NC004148_ 73_96_R	TGAACCAATTGCATTAGTCTC ACT	280
390	MNVL_NC004148_ 30_47_F	TGTCTATCTTCCTGACTC	157	MNVL_NC004148_ 73_91_R	TAATTGCATTAGTCTCACT	190
391	MNVL_NC004148_ 1907_1931_F	TAAGTAAGTTCAACCAAGCCTTTAG	5	MNVL_NC004148_ 1984_2006_R	TGTAACCAGCAGAATAGGCTT TG	332
392	MNVL_NC004148_ 1917_1931_F	TAACCAAGCCTTTAG	3	MNVL_NC004148_ 1984_1999_R	TACAGAATAGGCTTTG	191
393	MNVL_NC004148_ 2389_2414_F	TCAAGAGATCTTCAGTTTATGAGTAA	55	MNVL_NC004148_ 2468_2488_R	TGTTTATCCATGGTCCCACTC	354
394	MNVL_NC004148_ 2389_2408_F	TCAAGAGATCTTCAGTTTAT	54	MNVL_NC004148_ 2479_2498_R	TCTAATATTGTGTTTATCCA	261
395	MNVL_NC004148_ 2738_2756_F	TGAACATACCAATGCAGTT	104	MNVL_NC004148_ 2794_2810_R	TCAGGGGTCCTTCTATA	230
409	CRVCP_AY219836_ 2_22_F	TGACGTAGCCAAGGAGGCGTT	109	CRVCP_AY219836_ 46_70_R	TGAGGATGTGTCCAAGATGGC TGCG	287
410	CRVCP_AY219836_ 45_69_F	TCGCAGCCATCTTGGCAACATCCTC	95	CRVCP_AY219836_ 136_160_R	TGCGGGAAAGGCGGGAGTTGA AGAT	303
411	CRVCP_AY219836_ 475_499_F	TCCAGGTACTTCACCCCCAAACCTG	73	CRVCP_AY219836_ 577_601_R	TGCCGAGGCCTATGTGGTCGA CATT	301
412	CRVRP_AY219836_ 24_48_F	TCAACCACATAAGAGGTGGCTGTTC	52	CRVRP_AY219836_ 84_108_R	TAGGGAGATTGGGAGCTCCCG TATT	209
413	CRVRP_AY219836_ 440_464_F	TGGCTGAACTTTTGAAAGTGAGCGG	139	CRVRP_AY219836_ 495_517_R	TCGGGCCCACTATGACGTGTA CA	255

Primer Pair			Forward SEQ ID			Reverse SEQ ID
Number	Forward Primer Name	Forward Sequence	NO:	Reverse Primer Name	Reverse Sequence	NO:
414	CRVRP_AY219836_ 658_683_F	TGGGATGATCTACTGAGACTGTGTGA	141	CRVRP_AY219836_ 730_754_R	TGGTAATCAAAATACTGCGGG CCAA	325
415	CRVRP_AY219836_ 775_796_F	TGGTACTCCTCAACTGCTGTCC	146	CRVRP_AY219836_ 843_869_R	TCCGTGGATTGTTCTGTAGCA GTCTTC	244
990	HIV1_NC001802_ 7344_7364_F	TAGCAGGAAGCACTATGGGCG	25	HIV1_NC001802_ 7413_7439_R	TCAGCAAATTGTTCTGCTGCT GCACTA	224
991	HIV1_NC001802_ 7340_7361_F	TGAGCAGCAGGAAGCACTATGG	111	HIV1_NC001802_ 7413_7438_R	TAGCAAATTGCTTTGCTGTTG CACTA	205
992	HIV1_NC001802_ 4983_5014_F	TGTGAATATCAAGCAGGACATAACAA GGTAGG	163	HIV1_NC001802_ 5104_5127_R	TGGTCTTCTGGGGGCTTGTTCC ATC	330
993	HIV1_NC001802_ 1089_1117_F	TATAATCCACCTATCCCAGTAGGAGA AAT	40	HIV1_NC001802_ 1178_1204_R	TTTGGTCCTTGTCTTATGTCC AGAATG	370
994	HIV2_NC001722_ 8414_8439_F	TCGAAAAAACCTCCAGGCAAGAGTCAC	93	HIV2_NC001722_ 8476_8500_R	TCTAAACGCACATCCCCATGA ATTT	259
995	HIV2_NC001722_ 8425_8450_F	TCAGGCAAGAGTCACTGCTATCGAGA	65	HIV2_NC001722_ 8476_8502_R	TGTCTAAACCCACATCCCCAT GAATTT	341
996	HIV2_NC001722_ 5050_5075_F	TGCCAGGGAGTAGTAGAAGCAATGAA	121	HIV2_NC001722_ 5169_5196_R	TGTCATATCCCCTATTCCTCC CCTTCTT	337
997	HIV2_NC001722_ 5050_5075_F	TGCCAGGGAGTAGTAGAAGCAATGAA	121	HIV2_NC001722_ 5156_5187_R	TCCTATTCCTCCCCTTCTTTT AAAATTCATGC	245
998	HTLV1_NC001436_ 7221_7245_F	TGCCAATCACTCATACAACCCCCAA	118	HTLV1_NC001436_ 7330_7353_R	TGGTCTGGAAAAGACAGGGTT GGG	329
999	HTLV1_NC001436_ 7094_7118_F	TCAGAGCATCAGATCACCTGGGACC	63	HTLV1_NC001436_ 7153_7177_R	TGAGGGGAGTCGAGGGATAAG GAAC	289
1000	HTLV1_NC001436_ 7388_7410_F	GGAGGCTCCGTTGTCTGCATGTA	2	HTLV1_NC001436_ 7489_7516_R	TCGTTTGTAGCGAACATTGGT GAGGAAG	257
1001	HTLV1_NC001436_ 7818_7843_F	TACTCTCACACGGCCTCATACAGTAC	15	HTLV1_NC001436_ 7925_7947_R	TGGGGCTCATGGTCATTGTCA TC	322
1002	HTLV1_NC001436_ 7340_7361_F	TCTTTTCCAGACCCCGGACTCC	103	HTLV1_NC001436_ 7404_7428_R	TGGGAAAGCTGGTAGAGGTAC ATGC	316
1003	HTLV1_NC001436_ 7131_7156_F	TCGTTATCGGCTCAGCTCTACAGTTC	99	HTLV1_NC001436_ 7211_7233_R	TGAGTGATTGGCGGGGTAAGG AC	291
1004	HTLV2_NC001488_ 8180_8200_F	TAGAGGCGGATGACAATGGCG	21	HTLV2_NC001488_ 8254_8279_R	TACTTGGGATTGTTTGTGTGA GACGG	203
1005	HTLV2_NC001488_ 7757_7779_F	TCACCAAGGTGCCTCTAAAACGA	58	HTLV2_NC001488_ 7840_7861_R	TCATTGTGGTGGGTAGGTCGT C	233
1006	HTLV2_NC001488_ 2435_2456_F	TGGACCATCATTGGAAGGGACG	133	HTLV2_NC001488_ 2516_2540_R	TGGTGTTTGGAGTGGCTATTG GCAG	331
1007	HTLV2_NC001488_ 3592_3616_F	TTATGTCCTTGGGGTCACCTACTGG	172	HTLV2_NC001488_ 3680_3704_R	TCTTGCTTTGACATGTTGTGG TGGA	277
1008	HTLV2_NC001488_ 2880_2904_F	TCCATACTTCGCCTTCACCATTCCC	74	HTLV2_NC001488_ 2989_3013_R	TGTTGAGGACGGCTGCTAATT GTTG	348
1009	HTLV2_NC001488_ 1896_1919_F	TTTCGTTGTGGCAAGGTAGGACAC	184	HTLV2_NC001488_ 1993_2015_R	TTTGAGTTGTGGGCAGTCCCT TT	369
1010	HTLV2_NC001488_ 1198_1222_F	TCACCACGCAATGCTTCCCTATCTT	59	HTLV2_NC001488_ 1268_1291_R	TCCTGCTTGATGGCCTGTAAG TCT	247

3-continue	E 3-continued
------------	---------------

<b>D</b>		mer Pairs for Identification			110 1124000	D
Primer Pair			Forward SEQ ID			Reverse SEQ ID
Number	Forward Primer Name	Forward Sequence	NO :	Reverse Primer Name	Reverse Sequence	NO:
1011	HTLV2_NC001488_ 5735_5758_F	TTGGTCCATGACTCCGACCTTGAA	183	HTLV2_NC001488_ 5847_5870_R	TGAGGCTGGATCTATCCACGC AAA	288
1012	HTLV2_NC001488_ 5241_5265_F	TCCGATGCACATTCACGGTTGGTAT	80	HTLV2_NC001488_ 5332_5356_R	TAGTCGTTGGTCCGTTGTTAG GGAA	215
1013	HCV_NC001433_66_ 91_F	TCTAGCCATGGCGTTAGTATGAGTGT	101	HCV_NC001433_ 121_145_R	TGTTCCGCAGACTACTATGGC TCTC	372
1014	HCV_NC001433_66_ 91_F	TCTAGCCATGGCGTTACTATGAGTGT	101	HCV_NC001433_ 146_167_R	TGGCAATTCCGGTGTACTCAC C	311
1015	HCV_NC001433_66_ 91_F	TCTAGCCATGGCGTTAGTATGAGTGT	101	HCV_NC001433_ 128_153_R	TACTCACCGGTTCCGCAGACC ACTAT	197
1016	HCV_NC001433_51_ 74_F	TTCACGCAGAAAGCGTCTAGCCAT	175	HCV_NC001433_ 121_145_2_R	TGTTCCGCAGACCACTATGGC TCTC	346
1017	HCV_NC001433_51_ 73_F	TTCACGCAGAAAGCGTCTAGCCA	174	HCV_NC001433_ 146_167_R	TGGCAATTCCGGTGTACTCAC C	311
1018	HCV_NC001433_51_ 73_F	TCACGCAGAAAGCGTCTAGCCA	174	HCV_NC001433_ 128_153_R	TACTCACCGGTTCCGCAGACC ACTAT	197
1019	HCV_NC001433_62_ 85_F	AGCGTCTAGCCATGGCGTTAGTAT	1	HCV_NC001433_ 128_153_R	TACTCACCGGTTCCGCAGACC ACTAT	197
1020	HCV_NC001433_ 227_248_F	TCGCAAGACTGCTAGCCGAGTA	94	HCV_NC001433_ 277_298_R	TCGCAAGCACCCTATCAGGCA G	251
1021	HCV_NC001433_ 671_698_F	TAGGTCGCGTAATTTGGGTAAGGTCA TC	37	HCV_NC001433_ 720_742_R	TGAATGTACCCCATGAGGTCG GC	284
1022	HCV_NC001433_ 8598_8623_F	TCCTTCACGGAGGCTATGACTAGGTA	92	HCV_NC001433_ 8674_8700_R	TCGACACATTGGAGGAGCATG ATGTTA	249
1023	WN_NC001563_ 8365_8390_F	TCAGTGAATATCACTAGCCAGGTGCT	66	WN_NC001563_ 8434_8463_R	TGTTCCACTTCCCAAGTTTAC ATCTTCCTC	345
1024	WN_NC001563_ 8654_8678_F	TGCTCCTCTCAAAACCATGGGACAC	129	WN_NC001563_ 8749_8771_R	TCAGGAGCTTTCGTGTCCACC TT	227
1025	WN_NC001563_ 9026_9055_F	TCATCTACAACATGATGGGAAAGAGA GAGA	69	WN_NC001563_ 9101_9121_R	TAGCTCCGAGCCACATGAACC	206
1026	WN_NC001563_ 10135_10164_F	TGGATAGAGGAGAATGAATGGATGGA AGAC	135	WN_NC001563_ 10216_10237_R	TCAGGCTGCCACACCAGATGT C	229
1027	WN_NC001563_ 9016_9043_F	TGTCACACTTGCATTTACAACATGAT GG	155	WN_NC001563_ 9088_9112_R	TCCACATGAACCAAATGGCTC TGCT	235
1028	WN_NC001563_ 5696_5721_F	TCAAGATGGGGAATGAGATTGCCCTT	56	WN_NC001563_ 5758_5783_R	TACTCCGTCTCGTACGACTT TCTGTT	198
1029	WN_NC001563_ 5687_5714_F	TGCCTAGTGTGAAGATGGGGAATGAG AT	124	WN_NC001563_ 5796_5826_R	TGTTGTGATAACAAAGTCCCA ATCATCGTTC	349
1030	WN_NC001563_ 128_152_F	TCCGGGCTGTCAATATGCTAAAACG	83	WN_NC001563_ 187_212_R	TCGATCAGGCTCAACATAGCC CTCTT	250
1031	WN_NC001563_ 5994_6017_F	TCCGTCGCAAGTTGGTGATGAGTA	84	WN_NC001563_ 6079_6104_R	TGCATGTTGATGTTGTCCAGC ATGAT	300
1032	WN_NC001563_ 3527_3552_F	TGATTGACCCTTTTCAGTTGGGCCTT	115	WN_NC001563_ 3591_3613_R	TGCTGATCTTGGCTGTCCACC TC	307
1245	HBV_X51970_320_ 342_F	TCAACCTCCAATCACTCACCAAC	53	HBV_X51970_379_ 402_R	TATATGATAAAACGCCGCAGA CAC	216

	Pri	mer Pairs for Identification	of Adve	entitious Contamin	ant Viruses	
rimer Pair	Forward Primer		Forward SEQ ID	Reverse Primer		Reverse SEQ ID
lumber	Name	Forward Sequence	NO :	Name	Reverse Sequence	NO :
1246	HBV_X51970_317_ 341_F	TCCCCAATCTCCAATCACTCACCAA	76	HBV_X51970_379_ 406_R	TAGGAATATGATAAAACGCCG CAGACAC	208
1247	HBV_X51970_311_ 335_F	TCGCAGTCCCCAATCTCCAATCACT	96	HBV_X51970_382_ 410_R	TGAAGAGGAATATGATAAAAC GCCGCAGA	281
1248	HBV_X51970_1375_ 1399_F	TGGCTGCTAGGCTGTGCTGCCAACT	140	HBV_X51970_1412_ 1436_R	TACGGGACGTAAACAAAGGAC GTCC	195
1249	HBV_X51970_1777_ 1798_F	TACTAGGAGGCTGTAGGCATAAATTG GT	13	HBV_X51970_1868_ 1886_R	TCAGGCACAGCTTGGAGGC	228
1250	HBV_X51970_183_ 203_F	TACCCCTGCTCGTGTTACAGG	11	HBV_X51970_228_ 254_R	TGTCTAGACTCTGTGGTATTG TGAGGA	342
1251	HBV_X51970_180_ 202_F	TAGGACCCCTGCTCGTGTTACAG	26	HBV_X51970_262_ 289_R	TGCTCCCCCTAGAAAATTGAG AGAAGTC	306
1252	HBV_X51970_374_ 392_F	TGGATGTGTCTGCGGCGTT	136	HBV_X51970_423_ 450_R	TCCAGAAGAACCAACAAGAAG ATGAGGC	236
1253	HBV_X51970_368_ 388_F	TATCGCTGGATGTGTCTGCGG	46	HBV_X51970_424_ 453_R	TAATCCAGAAGAACCAACAAG AAGATGAGG	189
1254	HBV_X51970_312_ 335_F	TGCAGTCCCCAATCTCCAATCACT	117	HBV_X51970_376_ 398_R	TGATAAAACGCCGCAGACACA TC	292
2293	HTLV_TAX_GENE_ NC_001436_7107_ 7130_F	TCACCTGGGACCCCATCGATGGAC	60	HTLV_TAX_GENE_ NC_001436_7169_ 7192_R	TCTCTGGGTGGGGAAGGAGGG GAG	264
2294	HTLV_TAX_GENE_ NC_001436_7094_ 7118_F	TCAGAGCATCAGATCACCTGGGACC	63	HTLV_TAX_GENE_ NC_001436_7203_ 7226_R	TTGGCGGGGTGAGGACCTTGA GGG	365
2295	HTLV_TAX_GENE_ NC_001436_6989_ 7017_F	TACCCAGTCTACGTGTTTGGCGACTG TGT	9	HTLV_TAX_GENE_ NC_001436_7109_ 7129_R	TCCATCGATGGGGTCCCAGGT	240
2408	ARENAS_NC004296_ 474_494_F	TGGTGTGGTGAGAGTTTGGGA	149	ARENAS_NC004296_ 520_540_R	TGGCATGGTGCCAAACTGATT	313
2409	ARENAS_NC004296_ 474_494_2_F	TGGIGTGGTGAGAGTTTGGGA	145	ARENAS_NC004296_ 520_540_2_R	TGGCATIGTGCCAAACTGATT	314
2410	ARENAS_NC004296_ 931_953_F	TGTCTTTCAGGAGATGGATGGCC	160	ARENAS_NC004296_ 982_1002_R	TGTGTTTTCCCAAGCTCTTCC	344
2411	ARENAS_NC004296_ 931_953_2_F	TGTCTTTCAGGIGAIGGATGGCC	161	ARENAS_NC004296_ 982_1000_R	TGTTTTCCCAIGCCCTCCC	355
2412	ARENAS_NC004293_ 459_479_F	TGGTGTTGTGAGGGTCTGGGA	150	ARENAS_NC004293_ 505_527_R	TGCTGGCATGAACCAAACTGA TT	308
2413	ARENAS_NC004293_ 459_479_F	TGGTGTTGTGAGGGTCTGGGA	150	ARENAS_NC004293_ 511_534_R	TGGTCAAAGCTGGCATGAACC AAA	327
2414	ARENAS_NC004293_ 459_479_F	TGGTGTTGTGAGGGTCTGGGA	150	ARENAS_NC004293_ 511_534_2_R	TGGTCAIAGCIGGCATGAACC AAA	328
2415	ARENAS_NC004293_ 459_479_F	TGGTGTTGTGAGGGTCTGGGA	150	ARENAS_NC004293_ 511_534P_R	TGGTCAIAGCIGGCATGAACp CpAAA	328
2416	ARENAS_NC004293_ 915_937_F	TGTCTCTCTGGAGATGGGTGGCC	159	ARENAS_NC004293_ 975_994_R	TCAACACTGGTGTTGTCCCA	221
2417	ARENAS_NC004293_ 915_937_F	TGTCTCTCTGGAGATGGGTGGCC	159	ARENAS_NC004293_ 975_994P_R	TCAACACTGGTGTpTpGTCpC pCA	221

TABLE	3-cont	inued
IADUE	3-00110	Turea

Primer		mer Pairs for Identification	Forward			Perrover
Pair			SEQ ID			Reverse SEQ ID
Number	Forward Primer Name	Forward Sequence	NO:	Reverse Primer Name	Reverse Sequence	NO:
2418	ARENAS_NC004293_ 915_937_2_F	TGTCTCTCIGGIGATGGITGGCC	158	ARENAS_NC004293_ 975_994_R	TCAACACTGGTGTTGTCCCA	221
2419	ARENAS_NC004293_ 915_937_2_F	TGTCTCTCIGGIGATGGITGGCC	158	ARENAS_NC004293_ 975_994P_R	TCAACACTGGTGTpTpGTCpC pCA	221
2420	ARENAS_NC004293_ 915_937P_F	TGTCTCTCIGGIGATGGITpGGCpCp	158	ARENAS_NC004293_ 975_994P_R	TCAACACTGGTGTpTpGTCpC pCA	221
2423	POL_NC003461_ 1620_1649_F	TAGTCTCAAAGAGAAAGAGATCAAAC AAGA	38	POL_NC003461_ 1747_1772_R	TCAATCTCTCCCTTAACCATC CCATT	222
2424	POL_NC003461_ 2107_2138_F	TGTACAATCTATGTAGGAGATCCTTA CTGTCC	152	POL_NC003461_ 2215_2244_R	TGTCCATAATTTCTGGCAATA ACCTTCTAT	340
2425	POL_NC003461_ 2253_2279_F	TATCAGTGCAATCCATCTAGCAGCTG T	45	POL_NC003461_ 2302_2329_R	TGGCTTGATTGTCTCCTTGAA CCATTGC	315
2426	POL_NC003461_ 2253_2279_F	TATCAGTGCAATCCATCTAGCAGCTG T	45	POL_NC003461_ 2320_2349_R	TACTCTTGTTGTCACAGCTAT AGCTTGATT	200
2427	POL_NC003461_ 2253_2279_F	TATCAGTGCAATCCATCTAGCAGCTG T	45	POL_NC003461_ 2320_2352_R	TGGTACTCTTGATGTCACAGC TATAGCTTGATT	326
2428	POL_NC003461_ 2256_2279_F	TAGTGCAATCCATCTAGCAGCTGT	39	POL_NC003461_ 2302_2324_R	TGATTGTCGCCTTGAACCATT GC	293
2429	POL_NC003461_ 2463_2488_F	TGAGACCATCATAAGTAGCAAGATGT	110	POL_NC003461_ 2497_2523_R	TCTCCCATCATAGTATATCCT TTTGCT	263
2430	POL_NC003461_ 2935_2960_F	TATAGGGTCATGAATCAAGAACCCGG	41	POL_NC003461_ 2980_3004_R	TGCATGAATAAGGGTCTGAAG CCCA	299
2431	POL_NC003461_ 2935_2960_2_F	TATAGGGTCATGAATCAAGAACCIGG	42	POL_NC003461_ 2980_3004_R	TGCATGAATAAGGGTCTGAAG CCCA	299
2432	POL_NC003461_ 3644_3668_F	TTGGATCAGCCACTGATGAAAGATC	182	POL_NC003461_ 3760_3783_R	TGCTTCCATCCACGATATCTC ATC	309
2433	POL_NC003461_ 3759_3782_F	TGATGACATATCGTGGATGGAAGC	114	POL_NC003461_ 3880_3912_R	TGCACTAGAAAACTTCATCTG GGTTGCAGTATC	295
2434	POL_NC003461_ 3759_3782_2_F	TGATGAAATATCGTGGATGGAAGC	113	POL_NC003461_ 3880_3912_2_R	TGCACTAGAGAATTTCATCTG GGTTGCAGTATC	296
2435	RUBPOL_NC003443_ 1636_1661_F	TGTGCATCTTACTCACTAAAGGAGAA	164	RUBPOL_NC003443_ 1708_1734_R	TTCTGCAATTACTTGACACG ACCTCAT	358
2436	RUBPOL_NC003443_ 2067_2093_F	TCAGTTGATGGGCCTACCTCATTTCT T	68	RUBPOL_NC003443_ 2143_2167_R	TAGGGTCACTTGGAGGATTGA AAGG	210
2437	RUBPOL_NC003443_ 2133_2156_F	TGTAGGTGATCCCTTCAATCCTCC	154	RUBPOL_NC003443_ 2251_2275_R	TTGCAGAGATGGAGATCATTG TCCA	359
2438	RUBPOL_NC003443_ 2237_2261P_F	TATGTCAAAAGCTGTpGGACpAATpG AT	50	RUBPOL_NC003443_ 2318_2341P_R	ТТGCTTGGTTATpCpACpCpC ТGAACpCA	360
2439	RUBPOL_NC003443_ 2249_2276_F	TGTGGACAATCATCTCCATTGCTGCA AT	167	RUBPOL_NC003443_ 2301_2324_R	TGGACCATACAGGCAACCCGA CAA	310
2440	RUBPOL_NC003443_ 3689_3713_F	TTGCGTCAATGGCTTATATCAAAGG	181	RUBPOL_NC003443_ 3754_3778_R	TATCCCCGAAAGCCCAGATAT ATAC	217
2441	PNVL_U50363_ 2094_2120_F	TACAGATTTCAGCAAGTTCAATCAAG C	7	PNVL_U50363_ 2161_2181_R	TTGTGCACCATGCAGTTCATC	367
2442	PNVL_U50363_ 2110_2135_F	TTCAATCAAGCATTTCGGTATGAAAC	173	PNVL_U50363_ 2161_2181_R	TTGTGCACCATGCAGTTCATC	367

	_ Pri	mer Pairs for Identification	of Adve	ntitious Contamin	nant Viruses	
Primer Pair			Forward SEQ ID			Reverse SEQ ID
Number	Forward Primer Name	Forward Sequence	NO:	Reverse Primer Name	Reverse Sequence	NO:
2443	PNVL_U50363_ 2587_2612_F	TCACGAGATCTGCAGTTTATGAGTAA	62	PNVL_U50363_ 2677_2704_R	TGAAGTCATCCAGTATAGTGT TTATCCA	283
2444	PNVL_U50363_ 2581_2612_F	TATATATCACGAGATCTGCAGTTTAT GAGTAA	43	PNVL_U50363_ 2666_2686_R	TGTTTATCCACGGTCCCACTC	353
2445	PNVL_U50363_ 2583_2609_F	TATATCACGTGATCTGCAGTTTATGA G	44	PNVL_U50363_ 2677_2704_R	TGAAGTCATCCAGTATAGTGT TTATCCA	283
2446	PNVL_U50363_ 2660_2687_F	TCCTAAGAGTGGGACCATGGATAAAC AC	86	PNVL_U50363_ 2740_2769_R	TAATAGACTTTCCCCTCTATA TTCTAATTC	187
2447	PNVL_U50363_ 2660_2684_F	TCCTAAGAGTGGGACCATGGATAAA	85	PNVL_U50363_ 2740_2769_R	TAATAGACTTTCCCCTCTATA TTCTAATTC	187
2448	PNVL_U50363_ 2660_2684_F	TCCTAAGAGTGGGACCATGGATAAA	85	PNVL_U50363_ 2743_2776_R	TACTGCATAATAGGCTTTCCC CTCTAAATTCTAA	202
2449	PNVL_U50363_ 2663_2687_F	TAAGAGTGGGACCATGGATAAACAC	4	PNVL_U350363_ 2740_2769_2_R	TAATAGGCTTTCCCCTCTATA TTCTAATTC	188
2450	PNVL_U50363_ 1868_1894_F	TCAGTGTAGGTAGAATGTTTGCAATG C	67	PNVL_U50363_ 1918_1946_R	TCAGCTATCAATTTCTCTGCC AATATTTG	225
2451	PNVL_U50363_ 1874_1898_F	TAGGTAGAATGTTTGCAATGCAACC	36	PNVL_U350363_ 1918_1946_R	TCAGCTATCAATTTCTCTGCC AATATTTG	225
2452	PNVL_U50363_ 1872_1894_F	TGTAGGTAGAATGTTTGCAATGC	153	PNVL_U50363_ 1918_1946_2_R	TCAGCTATCATTTTCTCAGCC AAGATTTG	226
2453	MVL_AF266286_ 2002_2027_F	TGCCTTAATTGGAGATATGAGACCAT	127	MVL_AF266286_ 2035_2060_R	TAGATTTCATTTAGCCTCTGT GCAAA	204
2454	MVL_AF266286_ 2002_2027_2_F	TGCCTGAATTGGAGATATGAGACCAT	125	MVL_AF266286_ 2125_2144_R	TCGGGAGGGCAATGAGGGTC	254
2467	PNVL_U50363_ 2660_2684_1F	TCCTAAGAGTGGGACCATGGATAAA	85	PNVL_U50363_ 2743_2776_1R	TACTGCATAATAGACTTTCCC CTCTATATTCTAA	201
2533	PAV_IMP_ NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_ NC001526_6321_ 6355_R	TCTGCAACCATTTGCAAATAA TCTGGATATTTGCA	268
2534	PAV_IMP_ NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_ NC001526_6324_ 6355_R	TCTGCAACCATTTGCAAATAA TCTGGATATTT	267
2535	PAV_IMP_ NC001526_6222_ 6253_2_F	TAGGATGGCGATATGGTTGACACAGG CTTTGG	28	PAV_IMP_ NC001526_6324_ 6355_2_R	TCGGCAGCCATTTGCAAATAA TCAGGATATTT	253
2536	PAV_IMP_ NC001526_6212_ 6238_F	TACTGTTATTCAGGATGGTGATATGG T	19	PAV_IMP_ NC001526_6324_ 6355_R	TCTGCAACCATTTGCAAATAA TCTGGATATTT	267
2537	PAV_A9_NC001526_ 5632_5655_F	TTCAGATGTCTGTGTGGCGGCCTA	176	PAV_A9_NC001526_ 5691_5720_R	TACATATTCATCCGTGCTTAC AACCTTAGA	192
2538	PAV_A9_NC001526_ 5632_5655_F	TTCAGATGTCTGTGTGGGGGGCCTA	176	PAV_A9_NC001526_ 5691_5720_R	TACATATTCATCCGTGCTTAC AACCTTAGA	192
2539	PAV_A9_NC001526_ 2688_2715_F	TGGAAATCCTTTTTTCTCAAGGACGTG GT	131	PAV_A9_NC001526_ 2773_2802_R	TAGTATTTTGTCCTGCCACGC ATTTAAACG	211
2540	PAV_A9_NC001526_ 1972_2003_F	TAGATGATAGTGACATTGCATATAAA TATGCA	23	PAV_A9_NC001526_ 2085_2112_R	TTTCTGCTCGTTTATAATGTC TACACAT	368
2541	PAV_A7_NC001357_ 2011_2036_F	TATGGTGCAGTGGGCATTTGATAATG	49	PAV_A7_NC001357_ 2096_2121_R	TTGCTTTTTAAAAATGCAGCT GCATT	361

TABLE	3-continued
-------	-------------

Primer Pairs for Identification of Adventitious Contaminant Viruses						
Primer Pair			Forward SEQ ID			Reverse SEQ ID
Number	Forward Primer Name	Forward Sequence	NO:	Reverse Primer Name	Reverse Sequence	NO:
2542	PAV_A7_NC001357_ 2011_2036_F	TATGGTGCAGTGGGCATTTGATAATG	49	PAV_A7_NC001357_ 2108_2135_R	ТАТТТGCCTGCCААТТGCTTT ТТААААА	218
2543	PAV_A7_NC001357_ 4507_4530_F	TCCACCTGTGGTTATTGAACCTGT	72	PAV_A7_NC001357_ 4607_4629_R	TCAAACCCAGAGGTGCCTGTA AA	220
2544	PAV_A7_NC001357_ 748_767_F	TGACGAACCACAGCGTCACA	108	PAV_A7_NC001357_ 875_895_R	TGCACACAACGGACACACAAA	294
2545	PAV_A7_NC001357_ 947_966_F	TCGGGATGTAATGGCTGGTT	98	PAV_A7_NC001357_ 1034_1057_R	TACCATGTCCGAACCTGTATC TGT	193
2546	PAV_A10_ NC000904_875_ 903_F	TCAGGATGGTTTTTGGTAGAGGCTAT AGT	64	PAV_A10_NC000904_ 997_1027_R	TGCCTGTGCTTCCAAGGAATT GTGTGTAATA	302
2547	PAV_A10_ NC000904_1000_ 1027_F	TACACACAATTCCTTGGAAGCACAGG CA	6	PAV_A10_NC000904_ 1055_1079_R	TTAGGTCCTGCACAGCCGCAT AATG	356
2548	PAV_A10_ NC000904_1222_ 1246_F	TGAACTAACGGACAGTGGATATGGC	105	PAV_A10_NC000904_ 1275_1296_R	TCGCCATGTCTCTCTACCTGC G	252
2549	POLYOMA_ NC001669_2685_ 2713_F	TCAATGTATCTTATCATGTCTGGGTC CCC	57	POLYOMA_NC001669_ 2774_2796_R	_TGAGAGTGGATGGGCAGCCTA TG	286
2550	POLYOMA_ NC001669_2691_ 2719_F	TATCTTATCATGTCTGGGTCCCCAGG AAG	47	POLYOMA_NC001669_ 2774_2796_R	_TGAGAGTGGATGGGCAGCCTA TG	286
2551	POLYOMA_ NC001669_2698_ 2719_F	TCATGTCTGGGTCCCCTGGAAG	70	POLYOMA_NC001669_ 2774_2796_R	_TGAGAGTGGATGGGCAGCCTA TG	286
2552	POLYOMA_ NC001669_2726_ 2749_F	TGTGCCCTCAAAAACCCTGACCTC	166	POLYOMA_NC001669_ 2774_2796_R	_TGAGAGTGGATGGGCAGCCTA TG	286
2553	POLYOMA_ NC001669_2726_ 2749_2_F	TGTGCCCTCAAAAACCCTAACCTC	165	POLYOMA_NC001669_ 2774_2796_R	_TGAGAGTGGATGGGCAGCCTA TG	286
2554	POLYOMA_ NC001669_2767_ 2788_F	TGCCATTCATAGGCTGCCCATC	122	POLYOMA_NC001669_ 2899_2917_R	TCTGGAACCCAGCAGTGGA	275
2555	POLYOMA_ NC001669_4496_ 4520_F	TGAGGCCTATAGCAGCTATAGCCTC	112	POLYOMA_ NC001669_4583_ 4606_R	TAGCTGAAATTGCTGCTGGAG AGG	207
2556	POLYOMA_ NC001669_4533_ 4558_F	TCCCTCTACAGTAGCAACGGATGCAA	79	POLYOMA_ NC001669_4633_ 4660_R	TGGGGGACCTAATTGCTACTG TATCTGA	323
2557	POLYOMA_ NC001669_4537_ 4558_F	TCTACAGTAGCAACGGATGCAA	100	POLYOMA_ NC001669_4621_ 4641_R	TGTATCTGAAGCTGCTGCTGC	336
2558	POLYOMA_ NC001669_2978_ 3000_F	TAGGGGCCCAACCCCATTTTCAT	35	POLYOMA_ NC001669_3084_ 3110_R	TCCAGACCCTGCAAAAAATGA GAACAC	238
2559	POLYOMA_ NC001669_ 2967_2986_F	TCACCTTTGCAAAGGGGGCCC	61	POLYOMA_ NC001669_3087_ 3115_R	TGGGTCCCTGATCCAACTAGA AATGAAAA	324

rimer	:		Forward			Reverse
Pair	Forward Primer		SEQ ID	Reverse Primer		SEQ ID
lumber		Forward Sequence	NO:	Name	Reverse Sequence	NO :
2560	POLYOMA_ NC001669_3392_ 3412_F	TACGGTCACAGCTTCCCACAT	12	POLYOMA_ NC001669_3424_ 3446_R	TCTAAATGAGGACCTGACCTG TG	260
2561	POLYOMA_ NC001669_3780_ 3803_F	TCCTAGAAGTAGAGGCAGCATCCA	87	POLYOMA_ NC001669_3816_ 3840_R	TGCTCCAGGAGGTGCAAATCA AAGA	305
2643	HTLV1_NC001436_ 9387_7410_F	TGGAGGCTCCGTTGTCTGCATGTA	134	HTLV1_NC001436_ 7489_7516_R	TCGTTTGTAGGGAACATTGGT GAGGAAG	258
2670	PAV_IMP_MOD_ NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_MOD_ NC001526_6321_ 6355_R	TCTGCAACCATTTGCAAATAA TCTGGATATTTICA	269
2671	PAV_IMP_MOD_ NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_MOD_ NC001526_6321_ 6355_2_R	TCTGCAACCATTTGIAAATAA TCTGGATATTTICA	271
2672	PAV_IMP_MOD_ NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_MOD_ NC001526_6321_ 6355_3_R	TCTGCAACCATTTIIAAATAA TCTGGATATTTICA	273
2673	PAV_IMP_MOD_ NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_MOD_ NC001526_6324_ 6355P_R	ТСТGCААССАТрТрТрGAAAA ТААТСТGGATATTT	265
2674	PAV_IMP_MOD_ NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_MOD_ NC001526_6324_ 6355_R	TCTGCAACCATTTGAAAATAA TCTGGATATTT	266
2675	PAV_IMP_MOD_	TAGGATGGTGATATGGTTTGATA-	30	PAV_IMP_MOD_	TCTGCAACCATTTGIAAATAA	270
	NC001526_6222_ 6253_F	CAGG CTTTGG		NC001526_6324_ 6355_2_R	TCTGGATATTT	
2676	PAV_IMP_MOD_ NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_MOD_ NC001526_6324_ 6355_3_R	TCTGCAACCATTTIIAAATAA TCTGGATATTT	272
2677	PAV_IMP_MOD NC001526_6222_ 6253_F	TAGGATGGTGATATGGTTGATACAGG CTTTGG	30	PAV_IMP_MOD_ NC001526_6324_ 6355_4_R	TCTGCAACCATTTIIAIATAA TCTGGATATTT	274
2678	PAV_IMP_MOD_ NC001526_6324_ 6222_6253_2_F	TAGGATGGTGATATGGTTGATACAGG CTITGG	29	PAV_IMP_MOD_ NC001526_6324_ 6355_2_R	TCTGCAACCATTTGIAAATAA TCTGGATATTT	270
2679	PAV_IMP_MOD_ NC001526_6222_ 6253_3_F	TAGGATGGTGATATGGTTGATACAGG ITITGG	31	PAV_IMP_MOD_ NC001526_6324_ 6355_3_R	TCTGCAACCATTTIIAAATAA TCTGGATATTT	272
2680	PAV_IMP_MOD_ NC001526_6222_ 6253_4_F	TAGGATGGTGATATGGTTGATACIGG ITITGG	32	PAV_IMP_MOD_ NC001526_6324_ 6355_4_R	TCTGCAACCATTTIIAIATAA TCTGGATATTT	274
2681	PAV_IMP_MOD_ NC001526_6212_ 6238_F	TACTGTTATTCAGGATGGTGATATGG T	19	PAV_IMP_MOD_ NC001526_6324_ 6355_2_R	TCTGCAACCATTTGIAAATAA TCTGGATATTT	270
2682	PAV_IMP_MOD_ NC001526_6212_ 6238_F	TACTGTTATTCAGGATGGTGATATGG T	19	PAV_IMP_MOD_ NC001526_6324_ 6355_3_R	TCTGCAACCATTTIIAAATAA TCTGGATATTT	272
2683	PAV_IMP_MOD_ NC001526_6212_ 6238_F	TACTGTTATTCAGGATGGTGATATGG T	19	PAV_IMP_MOD_ NC001526_6324_ 6355_4_R	TCTGCAACCATTTIIAIATAA TCTGGATATTT	271

TABLE	3-continued

Primer Pair	:		Forward SEQ ID	l		Reverse SEQ ID
Number	Forward Primer Name	Forward Sequence	NO:	Reverse Primer Name	Reverse Sequence	NO:
2684	PAV_IMP_MOD_ NC001526_6212_ 6238_2_F	TACTGTTATICAGGATGGTGATATGG T	17	PAV_IMP_MOD_ NC001526_6324_ 6355_2_R	TCTGCAACCATTTGIAAATAA TCTGGATATTT	270
2685	PAV_IMP_MOD_ NC001526_6212_ 6238_3_F	TACTGTTATTCAGGATGGIGATATGG T	18	PAV_IMP_MOD_ NC001526_6324_ 6355_3_R	TCTGCAACCATTTIIAAATAA TCTGGATATTT	272
2686	PAV_IMP_MOD_ NC001526_6212_ 6238_4_F	TACTGTTATICAGGATGGIGATATGG T	16	PAV_IMP_MOD_ NC001526_6324_ 6355_4_R	TCTGCAACCATTTIIAIATAA TCTGGATATTT	274
2687	PAV_A9_MOD_ NC001526_5632_ 5655_F	TTCAGATGTCTGTGTGGGCIGCCTA	177	PAV_A9_MOD_ NC001526_5691_ 5720_R	TACATATTCATCCGTGCTTAC AACCTTAGA	192
2688	PAV_A9_MOD_ NC001526_5632_ 5655_2_F	TTCAGATGTCTITGTGGCIGCCTA	178	PAV_A9_MOD_ NC001526_5691_ 5720_R	TACATATTCATCCGTGCTTAC AACCTTAGA	192
2689	PAV_A9_MOD_ NC001526_2688_ 2715_F	TGGAAATCCTTTTTCTCAAGGACGTG GT	131	PAV_A9_MOD_ NC001526_2773_ 2802_R	TAGTATTTTGTCCTGCCACIC ATTTAAACG	212
2690	PAV_A9_MOD_ NC001526_2688_ 2715_F	TGGAAATCCTTTTTCTCAAGGACGTG GT	131	PAV_A9_MOD_ NC001526_2773_ 2802_2_R	TAGTATTTTGTCCTGCCAIIC ATTTAAACG	213
2691	PAV_A9_MOD_ NC001526_2688_ 2715_F	TGGAAATCCTTTTTCTCAAGGACGTG GT	131	PAV_A9_MOD_ NC001526_2773_ 2802_3_R	TAGTATTTTGTCCTGCCIIIC ATTTAAACG	214
2692	PAV_A9_MOD_ NC001526_1972 2003_F	TAGATGATAGTGAIATIGCATATIAA TATGCA	24	PAV_A9_MOD_ NC001526_2085_ 2112_R	TTTCTGCTCGTTTATAATGTC TACACAT	368
2693	PAV_A7_A10_ NC000904_1912_ 1937_F	TATGGTGCAGTGGGCATTTGATAATG	49	PAV_A7_A10_ NC000904_1997_ 2022_R	TTGCTTTTTAAAAATGCAGII GCATT	363
2694	PAV_A7_A10_ NC000904_1912_ 1937_F	TATGGTGCAGTGGGCATTTGATAATG	49	PAV_A7_A10_ NC000904_1997_ 2022_2_R	TTGCTTTTTAAAAATGCIIII GCATT	364
2695	PAV_A7_A10_ NC000904_1912_ 1937_2_F	TATGGTGCAGTGGGCATITGATAATG	48	PAV_A7_A10_ NC000904_1997_ 2022_R	TTGCTTTTTAAAAATGCAGII GCATT	362
2696	PAV_A7_A10_ NC000904_1912_ 1937_F	TATGGTGCAGTGGGCATTTGATAATG	49	PAV_A7_A10_ NC000904_2009_ 2036_R	TATTTGCCTGCIIATTGCTIT TTAAAAA	219
2807	PYV_NO_MAMMAL_ NC001663_3132_ 3159_F	TAGGGATTTTTGACCCATCTTTTTCT CA	34	PYV_NO_MAMMAL_ NC001663_3261_ 3281_R	TGGGCCTCTCTGCAAAGGAGA	319
2808	PYV_NO_MAMMAL_ NC001663_3267_ 3295_F	TTACAAAGAGGCCCAACGCCATTCTC ATC	171	PYV_NO_MAMMAL_ NC001663_3350_ 3386_R	TGAAAACACAAGATACTTTGG AACATACACAGGAGGT	279
2809	PYV_NO_MAMMAL_ NC001663_3578_ 3598_F	TCCTCCCACAGCAAACATGTG	88	PYV_NO_MAMMAL_ NC001663_3689_ 3716_R	TGTGTCTGTAAAGACTGAAGT TGTTGGA	343
2810	PYV_NO_MAMMAL_ NC001663_3561_ 3586_F	TGTAAATCTAGTGGCTCTCCTCCCAC	151	PYV_NO_MAMMAL_ NC001663_3686_ 3716_R	TGTAACTGTTAAAACTGAGGT TGTTGGAGTG	333

TABLE	3-continued	

Primer			Forward			Reverse
Pair			SEQ ID			SEQ II
Number	Forward Primer Name	Forward Sequence	NO:	Reverse Primer Name	Reverse Sequence	NO:
2811	PYV_NO_MAMMAL_ NC001663_3267_ 3294_F	TTGCAAAGAGGCCCAACCCCATTCTC AT	179	PYV_NO_MAMMAL_ NC001663_3351_ 3386_R	TGAGAACACAAGATACTITGG AAACTICACAGGIGG	285
2812	PYV_NO_MAMMAL_ NC001663_3267_ 3294_2_F	TTGCAAAGAGGCCCAACCCCATTITC AT	180	PYV_NO_MAMMAL_ NC001663_3351_ 3386_R	TGAGAACACAAGATACTITGG AAACTICACAGGIGG	285
2813	PYV_NO_MAMMAL_ NC001663_3267_ 3294_F	TTGCAAAGAGGCCCAACCCCATTCTC AT	179	PYV_NO_MAMMAL_ NC001663_3372_ 3404_R	TCCAGACCCAICCAAGAATGA GAACACAAGATA	237
2814	PYV_NO_MAMMAL_ NC001663_3267_ 3294_2_F	TTGCAAAGAGGCCCAACCCCATTITC AT	180	PYV_NO_MAMMAL_ NC001663_3372_ 3404_R	TCCAGACCCAICCAAGAATGA GAACACAAGATA	237
2864	AAV_NS1_ NC002077_1005_ 1028_F	TGCCACGACCGGCAAGACCAACAT	119	AAV_NS1_ NC002077_1126_ 1149_R	TGTCATCTTGCCCTCCTCCCA CCA	338
2865	AAV_NS1_ NC002077_1066_ 1094_F	TGCGTTAACTGGACCAATGAGAACTT TCC	128	AAV_NS1_ NC002077_1126_ 1149_2_R	TGTCATTTTGCCCTCCTCCCA CCA	339
2866	AAV_NS1_ NC002077_1020_ 1046_F	TACCAACATCGCGGAGGCTATIGCCC A	8	AAV_NS1_ NC002077_1126_ 1149_2_R	TGTCATTTTGCCCTCCTCCCA CCA	339
2867	AAV_NS1_ NC002077_1020_ 1046_F	TACCAACATCGCGGAGGCTATIGCCC A	8	AAV_NS1_ NC002077_1072_ 1099_R	TGAAGGGAAAGTTCTCATTGG TCCAGTT	282
2868	AAV_VP1_ NC002077_739_ 764_F	TGGGTCCTGCCCACCTACAACAACCA	144	AAV_VP1_ NC002077_832_ 859_R	TGTTGAAGTCAAAATACCCCC AGGGGGT	347
2869	AAV_VP1_ NC002077_743_ 764_F	TCCTGCCCACCTACAACAACCA	91	AAV_VP1_ NC002077_847_ 872_R	TGGCAGTGGAATCGGTTGAAG TCAAA	312
2870	AAV_VP1_ NC002077_648_ 665_F	TGGTGCCGATGGAGTGGG	148	AAV_VP1_ NC002077_680_ 703_R	TCCATTTGGAATCGCAATGCC AAT	242
2871	AAV_VP1_ NC002077_648_ 665_2_F	TGGTGCCGACGGAGTGGG	147	AAV_VP1_ NC002077_680_ 703_2_R	TCCATGGGGAATCGCAATGCC AAT	241
2872	AAV_VP1_ NC002077_831_ 851_F	TACCCCCTGGGGGTACTTTGA	10	AAV_VP1_ NC002077_886_ 910_R	TGTTGTTGATGAGTCTCTGCC AGTC	352
2873	AAV_VP1_ NC002077_838_ 867_F	TGGGGGTATTTTGACTTCAACCGATT CCAC	142	AAV_VP1_ NC002077_886_ 910_R	TGTTGTTGAPGAGTCTCTGCC AGTC	352
2874	AAV_VP1_ NC002077_838_ 867_2_F	TGGGGGTATTTTGACTTCAACIGITT CCAC	143	AAV_VP1_ NC002077_886_ 910_2_R	TGTTGTTGATGAGTC ICTGCC AGTC	351
2875	AAV_VP1_ NC002077_670_ 693_F	TCCTCGGGAAATTGGCATTGCGAT	89	AAV_VP1_ NC002077_748_ 772_R	TGTAGAGGTGGTTGTTGTAGG TGGG	335
2974	EBNA-2_NC007605- 36216_37679_780_ 803_F	TCCTCTCACTCATCAGAGCACCCC	90	EBNA-2_NC007605- 36216-37679_853_ 879_R	TGAGGGGGGATAATGGCATAGG AGGAAT	290

Primer Pair			Forward SEQ ID			Reverse SEQ ID
Number	Forward Primer Name	Forward Sequence	NO:	Reverse Primer Name	Reverse Sequence	NO:
2975	EBNA-2_NC007605- 36216-37679_2_ 29_F	TGCCTACATTCTATCTTGCGTTACAT GG	123	EBNA-2NC007605- 36216-37679_70_ 90_R	TACGGAGAGTGACGGGTTTCC	194
2976	EBNA-2_NC007605- 36216-37679_ 1195_1214_F	TACTCACCAACTCCTGGCCC	14	EBNA-2_NC007605- 36216-37679_ 1270_1292_R	TGGGACTCTGGTTCATGTATT GG	317
2977	EBNA-2_NC007605- 36216-37679_ 1319_1340_F	TCCCCGGTGATTGGTATCCTCC	78	EBNA-2_NC007605- 36216-37679_ 1390_1418_R	TCTTCATCAGAGCTAGGAGAT TCTGTTGT	276
2978	EBNA-3A_ NC007605-	TCCGGCCCTGGATGACAA	82	EBNA-3A_	TCTATGAGGACATTTTCCCAA	262
	NC007805- 79955-82877_24_ 41_F			NC007605-79955- 82877_94_119_R	тстсс	
2979	EBNA-3A_	TCGGCGCAAGTCCCAGAACC	97	EBNA-3A_	TACGGGGCCATGCCGTGTTG	196
	NC007605- 79955-82877_ 1318_1337_F			NC007605-79955- 82877_1387_ 1406_R		
2980	EBNA-3A_ NC007605-	TGGCCCCGTGTCTGGTAGC	137	EBNA-3A_	TCGGGTACTGGTGCACACGC	256
	79955-82877_ 1397_1415_F			NC007605-79955- 82877_1492_ 1511_R		
2981	EBNA-3A_ NC007605-79955- 82877_1617_1637_ F	TGTCCCGACTGTGGCACTTGA	156	EBNA-3A_ NC007605-79955- 82877_1669_ 1691_R	TGCATAGCAATCTCAGGAGGT GC	298
2982	EBNA-3B_ NC007605-83065- 85959_1950_1973_ F	TCCCCATCAGACACCTCAGGTGGA	77	EBNA-3B_ NC007605-83065- 85959_2002_ 2029_R	TGGGCTGATATGGAATGTGCC CTATCTG	320
2983	EBNA-3B_ NC007605-83065- 85959_891_ 911_F	TGAATGGTTCCGCCAGTGCAC	106	EBNA-3B_ NC007605-83065- 85959_946_968_R	TAAGGAACGGCGACAGGATGC GC	186
2984	EBNA-3B_ NC007605-83065- 85959_1046_1073_ F	TGCCAGATGATCCTATAATTGTTGAG GA	120	EBNA-3B_ NC007605-83065- 85959_1150_1171_ R	TGTACGGCAGTGTTTGCGGTA T	334
2985	EBNA-3B_ NC007605-83065- 85959_2531_2552_ F	TGACACAGGCACCCACGGAATA	107	EBNA-3B_ NC007605-83065- 85959_2593_ 2617_R	TCACTCGTTTAGACGGCGGAA TATC	223
2986	EBNA-3C_ NC007605-86083- 89135_1284_1307_ F	TGTGAAGCGCACAATTGTTAAGAC	162	EBNA-3C_ NC007605-86083- 89135_1330_ 1349_R	TCAGGGGTGCTTTGTGCTTC	231
2987	EBNA-3C_ NC007605-86083- 89135_2549_2573_ F	TCCCACATCTGCAATCGGAGACAGG	75	EBNA-3C_ NC007605-86083- 89135_2620_ 2642_R	TGGGGAGATGACCATGATGGT GC	321
2988	EBNA-3C_ NC007605-86083- 89135_1538_1562_ F	TGATTGATGTTGAAACCACCGAAGA	116	EBNA-3C_ NC007605-86083- 89135_1579_ 1600_R	TCCGATGTGGCTTATTTGGCT G	243

TABLE	3-continued
-------	-------------

TABLE	3-cont	inued
-------	--------	-------

	Pri	mer Pairs for Identification	of Adve	ntitious Contamin	ant Viruses	
Primer Pair Number	Forward Primer	Forward Sequence	Forward SEQ ID NO:	Reverse Primer Name	Reverse Sequence	Reverse SEQ ID NO:
2989	EBNA-1_NC007605- 95662-97587_ 1466_1493_F	TGGCTAGGTGTCACGTAGAAAGGACT AC	138	EBNA-1_NC007605- 95662-97587_ 1510_1533_R	TCCATATACGAACACACCGGC GAC	239
2990	EBNA-1_NC007605- 95662-97587 1420_1448_F	TCCAACCAGAAATTTGAGAACATTGC AGA	71	EBNA-1_NC007605- 95662-97587_ 1477_1499_R	TCCTCGGTAGTCCTTTCTACG TG	246
2991	LMP1_S75235-1- 1149_1023_1043_F	TCCGCCTTCGATGACACACGG	81	LMP1_S75235-1- 1149_1057_1076_ R	TGCAGCGTAGGAAGGTGTGG	297
2992	LMP1_S75235-1- 1149_600_620_F	TCTCCCGCACCCTCAACAAGC	102	LMP1_S75235-1- 1149_658_679_R	TCAGGTGGTGTCTGCCCTCGT T	232

**[0215]** Table 4 indicates the primer pair name virus identifier for the primer pairs disclosed herein.

# TABLE 4

Primer Pair Name Identifiers for Selected Viruses					
Virus Species	Virus Family	Primer Pair Name Virus Identifier	GenBank Accession Numbers		
Arenavirus	Arenaviridae	ARENAS	NC_004296		
Circovirus	Circoviridae	CRVCP	AY219836		
Hepatitis C virus	Flaviviridae	HCV	NC_001433		
West Nile virus	Flaviviridae	WN	NC_001563		
Hepatitis B virus	Hepadnaviridae	HBV	X51970		
Papillomavirus	Papillomaviridae	PAV_IMP	NC_001526		
Papillomavirus	Papillomaviridae	PAV_A9	NC_001526		
Papillomavirus	Papillomaviridae	PAV_A7	NC_001357		
Papillomavirus	Papillomaviridae	PAV_A10	NC_000904		
Respirovirus	Paramyxoviridae	RVL	X03614		
Pneumovirus	Paramyxoviridae	PVL	U50363		
Pneumovirus	Paramyxoviridae	PVNL	U50363		
Morbillivirus	Paramyxoviridae	MVL	AF266286		
Morbillivirus	Paramyxoviridae	MSVL	AF266286		
Metapneumovirus	Paramyxoviridae	MNVL	NC_004148		
Respirovirus	Paramyxoviridae	POL	NC_003461		
Rubulavirus	Paramyxoviridae	RUBPOL	NC_003443		
Human Immunodeficiency	Retroviridae	HIV1	NC_001802		
virus 1					
Human Immunodeficiency virus 2	Retroviridae	HIV-2	NC_001722		
Human T-lymphotropic virus 1	Retroviridae	HTLV1	NC_001436		
Human T-lymphotropic virus 2	Retroviridae	HTLV2	NC_001488		
Human T-cell Lymphotropic Virus	Retroviridae	HTLV_TAX_GENE	NC_001436		
Polyomavirus	Polyomaviridae	POLYOMA	NC_001669		
Non-mammalian	Polyomaviridae	PYV_NO_MAMMAL	NC_001663		
Polyomavirus	, 01111111111111		2.3_001005		
Dependovirus	Parvoviridae	AAV_NS1	NC 002077		
Dependovirus	Parvoviridae	AAV_VP1	NC_002077		
Dependovnus	1 ai voviniuae	AAY_Y11	110_002077		
#### Example 2

#### Sample Preparation and PCR

**[0216]** Samples were processed to obtain viral genomic material using a Qiagen QIAamp Virus BioRobot MDx Kit. Resulting genomic material was amplified using an Eppendorf thermal cycler and the amplicons were characterized on a Bruker Daltonics MicroTOF instrument. The resulting data was analyzed using GenX software (SAIC, San Diego, Calif. and Ibis, Carlsbad, Calif.).

**[0217]** All PCR reactions were assembled in 50  $\mu$ L reaction volumes in a 96-well microtiter plate format using a Packard MPII liquid handling robotic platform and M.J. Dyad thermocyclers (MJ research, Waltham, Mass.). The PCR reaction mixture consisted of 4 units of Amplitaq Gold, 1× buffer II (Applied Biosystems, Foster City, Calif.), 1.5 mM MgCl<sub>2</sub>, 0.4 M betaine, 800  $\mu$ M dNTP mixture and 250 nM of each primer. The following typical PCR conditions were used: 95° C. for 10 min followed by 8 cycles of 95° C. for 30 seconds, 48° C. for 30 seconds, and 72° C. 30 seconds with the 48° C. annealing temperature. increasing 0.9° C. with each of the eight cycles. The PCR was then continued for 37 additional cycles of 95° C. for 15 seconds, 56° C. for 20 seconds, and 72° C. 20 seconds.

#### Example 3

#### Solution Capture Purification of PCR Products for Mass Spectrometry with Ion Exchange Resin-Magnetic Beads

**[0218]** For solution capture of nucleic acids with ion exchange resin linked to magnetic beads,  $25 \ \mu$ l of a 2.5 mg/mL suspension of BioClone amine terminated superparamagnetic beads were added to 25 to 50  $\mu$ l of a PCR (or RT-PCR) reaction containing approximately 10 pM of a typical PCR amplification product. The above suspension was mixed for approximately 5 minutes by vortexing or pipetting, after which the liquid was removed after using a magnetic separator. The beads containing bound PCR amplification product were then washed three times with 50 mM ammonium bicarbonate/50% MeOH or 100 mM ammonium bicarbonate/50% MeOH. The bound PCR amplicon was eluted with a solution of 25 mM piperidine, 25 mM imidazole, 35% MeOH which included peptide calibration standards.

#### Example 4

#### Mass Spectrometry and Base Composition Analysis

**[0219]** The ESI-FTICR mass spectrometer is based on a Bruker Daltonics (Billerica, Mass.) Apex II 70e electrospray ionization Fourier transform ion cyclotron resonance mass spectrometer that employs an actively shielded 7 Tesla superconducting magnet. The active shielding constrains the majority of the fringing magnetic field from the superconducting magnet to a relatively small volume. Thus, components that might be adversely affected by stray magnetic fields, such as CRT monitors, robotic components, and other electronics, can operate in close proximity to the FTICR spectrometer. All aspects of pulse sequence control and data acquisition were performed on a 600 MHz Pentium II data station running Bruker's Xmass software under Windows NT 4.0 operating system. Sample aliquots, typically 15  $\mu$ l

were extracted directly from 96-well microtiter plates using a CTC HTS PAL autosampler (LEAP Technologies, Carrboro, N.C.) triggered by the FTICR data station. Samples were injected directly into a 10 µl sample loop integrated with a fluidics handling system that supplies the 100 µl/hr flow rate to the ESI source. Ions were formed via electrospray ionization in a modified Analytica (Branford, Conn.) source employing an off axis, grounded electrospray probe positioned approximately 1.5 cm from the metalized terminus of a glass desolvation capillary. The atmospheric pressure end of the glass capillary was biased at 6000 V relative to the ESI needle during data acquisition. A counter-current flow of dry N<sub>2</sub> was employed to assist in the desolvation process. Ions were accumulated in an external ion reservoir comprised of an rf-only hexapole, a skimmer cone, and an auxiliary gate electrode, prior to injection into the trapped ion cell where they were mass analyzed. Ionization duty cycles greater than 99% were achieved by simultaneously accumulating ions in the external ion reservoir during ion detection. Each detection event consisted of 1 M data points digitized over 2.3 s. To improve the signal-to-noise ratio (S/N), 32 scans were co-added for a total data acquisition time of 74 s.

**[0220]** The ESI-TOF mass spectrometer is based on a Bruker Daltonics MicroTOF<sup>TM</sup>. Ions from the ESI source undergo orthogonal ion extraction and are focused in a reflectron prior to detection. The TOF and FTICR are equipped with the same automated sample handling and fluidics described above. Ions are formed in the standard MicroTOF<sup>TM</sup> ESI source that is equipped with the same off-axis sprayer and glass capillary as the FTICR ESI source. Consequently, source conditions were the same as those described above. External ion accumulation was also employed to improve ionization duty cycle during data acquisition. Each detection event on the TOF was comprised of 75,000 data points digitized over 75  $\mu$ s.

[0221] The sample delivery scheme allows sample aliquots to be rapidly injected into the electrospray source at high flow rate and subsequently be electrosprayed at a much lower flow rate for improved ESI sensitivity. Prior to injecting a sample, a bolus of buffer was injected at a high flow rate to rinse the transfer line and spray needle to avoid sample contamination/carryover. Following the rinse step, the autosampler injected the next sample and the flow rate was switched to low flow. Following a brief equilibration delay, data acquisition commenced. As spectra were coadded, the autosampler continued rinsing the syringe and picking up buffer to rinse the injector and sample transfer line. In general, two syringe rinses and one injector rinse were required to minimize sample carryover. During a routine screening protocol a new sample mixture was injected every 106 seconds. More recently a fast wash station for the syringe needle has been implemented which, when combined with shorter acquisition times, facilitates the acquisition of mass spectra at a rate of just under one spectrum/minute.

**[0222]** Raw mass spectra were post-calibrated with an internal mass standard and deconvoluted to monoisotopic molecular masses. Unambiguous base compositions were

derived from the exact mass measurements of the complementary single-stranded oligonucleotides. Quantitative results are obtained by comparing the peak heights with an internal PCR calibration standard present in every PCR well at 500 molecules per well. Calibration methods are commonly owned and disclosed in U.S. Provisional Patent Application Ser. No. 60/545,425 which is incorporated herein by reference in entirety.

#### Example 5

#### De Novo Determination of Base Composition of Amplification Products Using Molecular Mass Modified Deoxynucleotide Triphosphates

[0223] Because the molecular masses of the four natural nucleobases have a relatively narrow molecular mass range (A=313.058, G=329.052, C=289.046, T=304.046-See Table 5), a persistent source of ambiguity in assignment of base composition can occur as follows: two nucleic acid strands having different base composition may have a difference of about 1 Da when the base composition difference between the two strands is  $G \cong A(-15.994)$  combined with C与T (+15.000). For example, one 99-mer nucleic acid strand having a base composition of A27G30C21T21 has a theoretical molecular mass of 30779.058 while another 99-mer nucleic acid strand having a base composition of  $A_{26}G_{31}C_{22}T_{20}$  has a theoretical molecular mass of 30780.052. A 1 Da difference in molecular mass may be within the experimental error of a molecular mass measurement and thus, the relatively narrow molecular mass range of the four natural nucleobases imposes an uncertainty factor.

**[0224]** The present invention provides for a means for removing this theoretical 1 Da uncertainty factor through amplification of a nucleic acid with one mass-tagged nucleobase and three natural nucleobases. The term "nucleobase" as used herein is synonymous with other terms in use in the art including "nucleotide,""deoxynucleotide,""nucleotide residue,""nucleotide triphosphate (NTP)," or deoxynucleotide triphosphate (dNTP).

[0225] Addition of significant mass to one of the 4 nucleobases (dNTPs) in an amplification reaction, or in the primers themselves, will result in a significant difference in mass of the resulting amplification product (significantly greater than 1 Da) arising from ambiguities arising from the  $G \leftrightarrows A$ combined with  $C \hookrightarrow T$  event (Table 5). Thus, the same the  $G \Leftrightarrow A$  (-15.994) event combined with 5-Iodo-C $\Leftrightarrow T$ (-110.900) event would result in a molecular mass difference of 126.894. If the molecular mass of the base composition  $A_{27}G_{30}$  5-Iodo- $C_{21}T_{21}$  (33422.958) is compared with  $A_{26}G_{31}\text{5-Iodo-}C_{22}T_{20}\text{, (33549.852)}$  the theoretical molecular mass difference is +126.894. The experimental error of a molecular mass measurement is not significant with regard to this molecular mass difference. Furthermore, the only base composition consistent with a measured molecular mass of the 99-mer nucleic acid is  $A_{27}G_{30}5$ -Iodo- $C_{21}T_{21}$ . In contrast, the analogous amplification without the mass tag has 18 possible base compositions.

TABLE 5

Molecular Masses of Natural Nucleobases and the Mass-Modified
Nucleobase 5-Iodo-C and Molecular Mass Differences
Resulting from Transitions

Nucleobase	Molecular Mass	Transition	Molecular Mass
A	313.058	A>T	-9.012
А	313.058	A>C	-24.012
А	313.058	A>5-Iodo-C	101.888
А	313.058	A>G	15.994
Т	304.046	T>A	9.012
Т	304.046	T>C	-15.000
Т	304.046	T>5-Iodo-C	110.900
Т	304.046	T>G	25.006
С	289.046	C>A	24.012
С	289.046	C>T	15.000
С	289.046	C>G	40.006
5-Iodo-C	414.946	5-Iodo-C>A	-101.888
5-Iodo-C	414.946	5-Iodo-C>T	-110.900
5-Iodo-C	414.946	5-Iodo-C>G	-85.894
G	329.052	G>A	-15.994
G	329.052	G>T	-25.006
G	329.052	G>C	-40.006
G	329.052	G>5-Iodo-C	85.894

**[0226]** Mass spectra of bioagent-identifying amplicons were analyzed independently using a maximum-likelihood processor, such as is widely used in radar signal processing. This processor, referred to as GenX, first makes maximum likelihood estimates of the input to the mass spectrometer for each primer by running matched filters for each base composition aggregate on the input data. This includes the GenX response to a calibrant for each primer.

[0227] The algorithm emphasizes performance predictions culminating in probability-of-detection versus probability-of-false-alarm plots for conditions involving complex backgrounds of naturally occurring organisms and environmental contaminants. Matched filters consist of a priori expectations of signal values given the set of primers used for each of the bioagents. A genomic sequence database is used to define the mass base count matched filters. The database contains the sequences of known bacterial bioagents and includes threat organisms as well as benign background organisms. The latter is used to estimate and subtract the spectral signature produced by the background organisms. A maximum likelihood detection of known background organisms is implemented using matched filters and a running-sum estimate of the noise covariance. Background signal strengths are estimated and used along with the matched filters to form signatures which are then subtracted. The maximum likelihood process is applied to this "cleaned up" data in a similar manner employing matched filters for the organisms and a running-sum estimate of the noisecovariance for the cleaned up data.

**[0228]** The amplitudes of all base compositions of bioagent-identifying amplicons for each primer are calibrated and a final maximum likelihood amplitude estimate per organism is made based upon the multiple single primer estimates. Models of all system noise are factored into this two-stage maximum likelihood calculation. The processor reports the number of molecules of each base composition contained in the spectra. The quantity of amplification product corresponding to the appropriate primer set is reported as well as the quantities of primers remaining upon completion of the amplification reaction. **[0229]** Base count blurring can be carried out as follows. "Electronic PCR" can be conducted on nucleotide sequences of the desired bioagents to obtain the different expected base counts that could be obtained for each primer pair. See for example, ncbi.nlm.nih.gov/sutils/e-pcr/; Schuler, Genome Res. 7:541-50, 1997. In one illustrative embodiment, one or more spreadsheets, such as Microsoft Excel workbooks contain a plurality of worksheets. First in this example, there is a worksheet with a name similar to the workbook name; this worksheet contains the raw electronic PCR data. Second, there is a worksheet named "filtered bioagents base count" that contains bioagent name and base count; there is a separate record for each strain after removing sequences that are not identified with a genus and species and removing all sequences for bioagents with less than 10 strains. Third, there is a worksheet, "Sheet1" that contains the frequency of substitutions, insertions, or deletions for this primer pair. This data is generated by first creating a pivot table from the data in the "filtered bioagents base count" worksheet and then executing an Excel VBA macro. The macro creates a table of differences in base counts for bioagents of the same species, but different strains. One of ordinary skill in the art may understand additional pathways for obtaining similar table differences without undo experimentation.

**[0230]** Application of an exemplary script, involves the user defining a threshold that specifies the fraction of the strains that are represented by the reference set of base counts for each bioagent. The reference set of base counts as are needed to meet or exceed the threshold. The set of reference base counts is defined by taking the most abundant strain's base type composition and adding it to the reference set and then the next most abundant strain's base type composition is added until the threshold is met or exceeded. The current set of data was obtained using a threshold of 55%, which was obtained empirically.

[0231] For each base count not included in the reference base count set for that bioagent, the script then proceeds to determine the manner in which the current base count differs from each of the base counts in the reference set. This difference may be represented as a combination of substitutions, Si=Xi, and insertions, Ii=Yi, or deletions, Di=Zi. If there is more than one reference base count, then the reported difference is chosen using rules that aim to minimize the number of changes and, in instances with the same number of changes, minimize the number of insertions or deletions. Therefore, the primary rule is to identify the difference with the minimum sum (Xi+Yi) or (Xi+Zi), e.g., one insertion rather than two substitutions. If there are two or more differences with the minimum sum, then the one that will be reported is the one that contains the most substitutions.

**[0232]** Differences between a base count and a reference composition are categorized as one, two, or more substitutions, one, two, or more insertions, one, two, or more deletions, and combinations of substitutions and insertions or deletions. The different classes of nucleobase changes and their probabilities of occurrence have been delineated in U.S. Patent Application Publication No. 2004209260 (U.S. application Ser. No. 10/418,514) which is incorporated herein by reference in entirety.

#### Example 6

#### Codon Base Composition Analysis—Assay Development

[0233] The information obtained by the codon analysis method of the present invention is base composition. While base composition is not as information-rich as sequence, it can have the same practical utility in many situations. The genetic code uses all 64 possible permutations of four different nucleotides in a sequence of three, where each amino acid can be assigned to as few as one and as many as six codons. Since base composition analysis can only identify unique combinations, without determining the order, one might think that it would not be useful in genetic analysis. However, many problems of genetic analysis start with information that constrains the problem. For example, if there is prior knowledge of the biological bounds of a particular genetic analysis, the base composition may provide all the necessary and useful information. If one starts with prior knowledge of the starting sequence, and is interested in identifying variants from it, the utility of base composition depends upon the codons used an the amino acids of interest.

[0234] Analysis of the genetic code reveals three situations, illustrated in Tables 6A-C. In Table 6A, where the leucine codon CTA is comprised of three different nucleotides, each of the nine possible single mutations are always identifiable using base composition alone, and result in either a "silent" mutation, where the amino acid is not changed, or an unambiguous change to another specific amino acid. Irregardless, the resulting encoded amino acid is known, which is equivalent to the information obtained from sequencing. In Table 6B, where two of the three nucleotides of the original codon are the same, there is a loss of information from a base composition measurement compared to sequencing. In this case, three of the nine possible single mutations produce unambiguous amino acid choices, while the other six each produce two indistinguishable options. For example, if starting with the phenylalanine codon TTC, then either one of the two Ts could change to A, and base composition analysis could not distinguish a first position change from a second position change. A first position change of T to A would encode an isoleucine and a second position change of T to A would encode a tyrosine. However no other options are possible and the value of the information would depend upon whether distinguishing an encoded isoleucine from a tyrosine was biologically important. In Table 6C, all three positions have the same nucleotide, and therefore the ambiguity in amino acid identity is increased to three possibilities. Out of 64 codon choices, 20 have three unique nucleotides (as in Table 6A), 40 have two of the same and one different nucleotide (as in Table 6B) and 4 have the same nucleotide in all three positions (as in Table 6C).

TABLE 6A

Wild Type Codon with Three Unique Nucleobases						
Codon Description	Codon(s)	Codon Base Composition	Amino Acid Coded			
WILD TYPE CODON Single Mutation Single Mutation Single Mutation Single Mutation	CTA ATA GTA TTA CAA	A1C1T1 A2T1 A1G1T1 A1T2 A1C2	Leu Ile Val Leu Gln			

Wild Type	Codon with T	hree Unique N	ucleobases
Codon		Codon Base	
Description	$\operatorname{Codon}(s)$	Composition	Amino Acid Coded
Single Mutation	CGA	A1G1C1	Arg
Single Mutation	CCA	A1C2	Pro
Single Mutation	CTG	G1C1T1	Leu
Single Mutation	CTC	C2T1	Leu
Single Mutation	CTT	C1T2	Leu

#### [0235]

37711

#### TABLE 6B

**.**...

**NT 1 1** 

Wild Type	Codon with 1	wo Unique Nu	cleobases
Codon Description	Codon(s)	Codon Base Composition	Amino Acid Coded
WILD TYPE CODON Single Mutations	TTC ATC, TAC	C1T2 A1C1T1	Phe Ile, Tyr
Single Mutations	GTC, TGC	G1C1T1	Val, Cys
Single Mutations	CTC, TCC	C2T1	Leu, Ser
Single Mutation	TTA	A1T2	Leu
Single Mutation	TTG	G1T2	Leu
Single Mutation	TTT	Т3	Phe

1.1 m

[0236]

TABLE 6C

Wild Type Co	don Having Three of	of the Same Nu	icleobase
Codon	Codon(s)	Codon Base	Amino
Description		Composition	Acid Coded
WILD TYPE CODON	TTT	T3	Phe
Single Mutations	ATT, TAT, TTA	A1T2	Ile, Tyr, Leu
Single Mutations	GTT, TGT, TTG	G1T2	Val, Cys, Leu
Single Mutations	CTT, TCT, TTC	C1T2	Leu, Ser, Phe

#### Example 7

#### Testing of HIV Strains for Drug Resistance by Codon Analysis

[0237] Upon determination of an HIV positive test, drug resistance testing can be done by conventional sequencing methods to establish the nucleotide sequence of the major HIV strain infecting the patient followed by analysis of the codons most important in mediating drug resistance. The patient would then be monitored while on antiretroviral therapy by codon analysis methods according to the present invention for the appearance of emerging viral mutations. The advantages of monitoring by these methods for rapid codon analysis are: (i) it is much more sensitive to identifying low abundance mutations in a population, (ii) it can be done on a much lower viral titer and (iii) it is less expensive than sequencing. Based upon previous data on the sensitivity and dynamic range of the mass spectrometer, the ability is anticipated to identify a low-abundance mutation present in as little as 0.1% of the viral population in a 10 to 100-fold lower titer of virus than can be analyzed by sequencing.

[0238] According to the 2005 update of drug resistance mutations in HIV there are 28 major and minor mutations in HIV reverse transcriptase and 21 major and minor mutations in HIV protease (Johnson, V. A. et al. 2005, Top. HV Med. 13, 51-57). A smaller subset of these mutations is currently used in making decisions regarding the drug regimen of choice for a particular virus strain. All 49 positions where mutations occur which are related to antiretroviral therapy were examined both with respect to the amino acid changes that mediate resistance and sequences present in isolates of the virus which included HIV-1 subtype B sequences from persons with well-characterized histories of antiretroviral treatment (Rhee, S. Y. et al. 2005, J. Infect. Dis. 192, 456-65). Sequences were obtained from the Stanford HIV Reverse Transcriptase and Protease Sequence Database and included sequences from published studies and previously unpublished sequences generated at Stanford University from patients living in northern California (GenBank accession numbers AY796421-AY798497 and AY800656-AY802758). For the present example, sequences were aligned using an alignment editor and the relevant codons of the reverse transcriptase gene in 2,102 sequences were analyzed. An example of this analysis is illustrated in Table 7. Out of 2,102 HIV sequences, the majority have a wild type codon encoding methionine in position 41 of reverse transcriptase. Deviations from wild type observed were all changes that encode leucine, which is a thymidine nucleotide-associated mutation (TAM) associated with drug resistance.

TABLE 7

	e					
Codon	Total	Fraction	n Amino Acid	Codon Base Composition	Phenotype	Distinguishable from Wild Type
ATG	1224	0.59	N (1) (1	L C T		
AIG	1224	0.58	Met[M]	$A_1G_1T_1$	Wild type	_
TTG	429	0.58	Met[M] Leu[L]	$A_1G_1T_1$ $G_1T_2$	Wild type TAM	yes
					21	

	TABLE 7-continued						
Mutations of Codon 41 in HIV Reverse Transcriptase							
Codon	Total	Fraction	Amino Acid	Codon Base Composition	Phenotype	Distinguishable from Wild Type	
CTA CTT	7		Leu[L] Leu[L]	$\begin{array}{c} A_1 C_1 T_1 \\ C_1 T_2 \end{array}$	TAM TAM	yes yes	
Subtotal Mixtures	1961 141	0.93 0.07					
Total	2102						

**[0239]** From codon analysis of this dataset, several observations can be made. First, the wild type codon is the most frequently observed (in this case 58%). Second, the next most abundant codons are single point mutations from wild type, and encode a documented mutation which causes drug-resistance. Third, a significant fraction of the sequences (in this case 7%) contain ambiguities with respect to the correct nucleotide in the sequencing reaction (in Tables 6B and 6C, all sequences with more than one codon per base composition are grouped together). These may be derived from the presence of a mixed population of virus sequences in the original patient samples. Codon 41 of reverse transcriptase exemplifies the situation in Table 6A where the codon is comprised of three different nucleotides

ATG, and all mutations are unambiguously identifiable with respect to the amino acid that they encode.

**[0240]** The reverse transcriptase codon leucine 210 (Table 8) provides an example of the wild type codon base composition situation illustrated in Table 6B. Leucine is the most extreme example of degeneracy in the genetic code, and six different codons encode leucine. Nevertheless, mutation of the wild type codon TTG to TGG (tryptophan), which causes drug resistance, is unambiguously distinguishable from all six wild type leucine codons present, and no other codon in the dataset contains the composition base composition  $G_2T$ . Therefore, base composition analysis may be sufficient to distinguish a drug-resistant strain from the wild type strain, even in the absence of prior knowledge of the particular codon used to encode leucine.

TABLE 8

Codon	Total	Fraction	n Amino Acid	Codon Base Composition	Phenotype	Distinguishable from Major or Minor Wild Typ
TTG	1223	0.58	Leu[L]	G1T2	Major wild type	_
TGG	497	0.24	Trp[W]	G2T1	TAM	yes
TTA	168	0.08	Leu[L]	A1T2	Minor Wild type	—
CTG	17	0.01	Leu[L]	G1C1T1	Minor Wild type	—
TTC	9	0.00	Phe[F]	C1T2		
TCA	5	0.00	Ser[S]	A1C1T1		
CTC	4	0.00	Leu[L]	C2T1	Minor Wild type	—
CTT	2	0.00	Leu[L]	C1T2	Minor Wild type	_
TCG	2	0.00	Ser[S]	G1C1T1		
СТА	1	0.00	Leu[L]	A1C1T1	Minor Wild type	_
GAG	1	0.00	Glu[E]	A1G2		
ATG	2	0.00	Met[M]	A1G1T1		
Subtotal	1931	0.92				
Mixtures	170	0.08				
Total	2101					

**[0241]** Not all mutations are unambiguously identifiable by base composition in the absence of knowledge of the starting sequence. For example, reverse transcriptase lysine 65 can mutate to arginine, which mediates drug resistance. Lysine is an example of the situation in **FIG. 2C**, where one of the two codons that encode lysine have all three positions comprised of the same nucleotide (AAA), and the second lysine codon has two of the same and one different nucleotide (AAG). A mutation that encodes arginine AGA is not distinguishable from the minor wild type codon, which does not allow unambiguous assignment of drug resistance using base composition. However, if the starting virus wild type sequence is known (either AAA or AAG), then any mutation to arginine would be distinguishable.

### TABLE 9

	Mutations of Codon 65 in HIV Reverse Transcriptase					
Codon	Total	Fraction	1 Amino Acid	Codon Base Composition	Phenotype	Distinguishable from Major or Minor Wild Type
AAA	1782	0.85	Lys[L]	A3	Wild type	
AAG	196	0.09	Lys[L]	A2G1	Minor	
					wild type	
AGA	46	0.02	Arg[R]	A2G1	NNRTI*	Not
						distinguishable
						from minor wild
ACC	1	0.00	4 [D]	1101	NUTIN	type
AGG	1	0.00	Arg[R]	A1G1	NNRTI*	Yes
subtotal	2025					
Mixtures	77	0.04				
		•				
total	2102					

\*NNRTI = nucleoside and nucleotide reverse transcriptase inhibitor resistant

**[0242]** Thus, for most of the 49 mutations cited in the 2005 update of drug resistance mutation in HIV (Johnson, V. A. et al. 2005, Top. HIV Med. 13, 51-57), base composition analysis is sufficient to determine whether a drug-resistant mutation is present, simply based upon the measured composition, the known biological constraints of the HIV virus and the rules of the genetic code. For a minority of the mutations, prior knowledge of the sequence of the particular virus strain is needed to distinguish mutations that mediate drug-resistance mutations from minor wild type codons variations.

**[0243]** The advantages of monitoring codon base compositions include: sensitivity to identification of low abundance mutations in a population, it can be done on a much lower viral titer and is less expensive than sequencing. Based upon previous data on the sensitivity and dynamic range of electrospray mass spectrometry, we would anticipate the ability to identify a low-abundance mutation present in as little as 0.1% of the viral population in a 10 to 100-fold lower titer of virus than can be analyzed by sequencing. If **[0245]** While in accordance with the patent statutes, description of the various embodiments and examples have been provided, the scope of the invention is not to be limited thereto or thereby. Modifications and alterations of the present invention will be apparent to those skilled in the art without departing from the scope and spirit of the present invention.

**[0246]** Therefore, it will be appreciated that the scope of this invention is to be defined by the appended claims, rather than by the specific examples which have been presented by way of example.

**[0247]** Each reference (including, but not limited to, journal articles, U.S. and non-U.S. patents, patent application publications, international patent application publications, gene bank accession numbers, internet web sites, and the like) cited in the present application is incorporated herein by reference in its entirety.

a drug-resistant mutation should arise as a small fraction of the population in a patient that currently has a low viral titer, it could be valuable to change drugs sooner rather than later to keep the viral titer low. Lower virus titers mean there are fewer viruses available to mutate, and drug resistance would be suppressed.

**[0244]** The present invention includes any combination of the various species and subgeneric groupings falling within the generic disclosure. This invention therefore includes the generic description of the invention with a proviso or negative limitation removing any subject matter from the genus, regardless of whether or not the excised material is specifically recited herein.

#### SEQUENCE LISTING

<160> NUMBER OF SEQ ID NOS: 372 <210> SEQ ID NO 1 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 1 agcgtctagc catggcgtta gtat 24 <210> SEQ ID NO 2 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 2 ggaggeteeg ttgtetgeat gta 23 <210> SEQ ID NO 3 <211> LENGTH: 15 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 3 taaccaagcc tttag 15 <210> SEQ ID NO 4 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 4 25 taagagtggg accatggata aacac <210> SEQ ID NO 5 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 5 taagtaagtt caaccaagcc tttag 25 <210> SEQ ID NO 6 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 6 28 tacacacaat tccttggaag cacaggca

-00		

<210> SEQ ID NO 7 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 7 tacagatttc agcaagttca atcaagc 27 <210> SEQ ID NO 8 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (22)...(22) <223> OTHER INFORMATION: I <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 8 taccaacatc gcggaggcta ttgccca 27 <210> SEQ ID NO 9 <211> LENGTH: 29 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 9 tacccagtct acgtgtttgg cgactgtgt 29 <210> SEQ ID NO 10 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 10 taccccctgg gggtactttg a 21 <210> SEQ ID NO 11 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 11 tacccctgct cgtgttacag g 21 <210> SEQ ID NO 12 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 12 21 tacggtcaca gcttcccaca t

<210> SEQ ID NO 13 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 13 tactaggagg ctgtaggcat aaattggt 28 <210> SEQ ID NO 14 <211> LENGTH: 20 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 14 tactcaccaa ctcctggccc 20 <210> SEQ ID NO 15 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 15 tactctcaca cggcctcata cagtac 26 <210> SEQ ID NO 16 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (10)...(10) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (19)...(19) <223> OTHER INFORMATION: I <400> SEQUENCE: 16 tactgttatt caggatggtg atatggt 27 <210> SEQ ID NO 17 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (10)...(10) <223> OTHER INFORMATION: I <400> SEQUENCE: 17 tactgttatt caggatggtg atatggt 27 <210> SEQ ID NO 18 <211> LENGTH: 27 <212> TYPE: DNA

-c	$\cap$	nt		n	11	0	а
$-\circ$	<u> </u>		_		u	c	u
<u> </u>	۰.		-	**	u	~	-

<213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (19)...(19) <223> OTHER INFORMATION: I <400> SEQUENCE: 18 tactgttatt caggatggtg atatggt 27 <210> SEQ ID NO 19 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 19 tactgttatt caggatggtg atatggt 27 <210> SEQ ID NO 20 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 20 tagaatgttt gctatgcaac c 21 <210> SEQ ID NO 21 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 21 tagaggcgga tgacaatggc g 21 <210> SEQ ID NO 22 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 22 tagatccaca agctttaggg tctg 24 <210> SEQ ID NO 23 <211> LENGTH: 32 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 23 tagatgatag tgacattgca tataaatatg ca 32 <210> SEQ ID NO 24 <211> LENGTH: 32 <212> TYPE: DNA

-co	ntinued
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<220> FEATURE:	
<221> NAME/KEY: modified_base	
<222> LOCATION: (14)(14)	
<223> OTHER INFORMATION: I	
<220> FEATURE:	
<221> NAME/KEY: modified_base	
<222> LOCATION: (17)(17)	
<223> OTHER INFORMATION: I	
<220> FEATURE:	
<221> NAME/KEY: modified_base	
<222> LOCATION: (24)(24) <223> OTHER INFORMATION: I	
<400> SEQUENCE: 24	
tagatgatag tgacattgca tataaatatg ca	32
<210> SEQ ID NO 25	
<211> LENGTH: 21	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 25	
tagcaggaag cactatgggc g	21
<210> SEQ ID NO 26	
<211> LENGTH: 23	
<211> HINGIN: 23 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 26	
taggacccct gctcgtgtta cag	23
<210> SEQ ID NO 27	
<211> LENGTH: 19	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 27	
taggaccgtg gataaacac	19
<210> SEQ ID NO 28	
<211> LENGTH: 32	
<211> LENGIN: 52 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 28	
taggatggcg atatggttga cacaggcttt gg	32
<210> SEQ ID NO 29	
<211> LENGTH: 32	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<220> FEATURE:	

	-continued
<pre>&lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (29)(29) &lt;223&gt; OTHER INFORMATION: I</pre>	
<400> SEQUENCE: 29	
taggatggtg atatggttga tacaggcttt gg	32
<pre>&lt;210&gt; SEQ ID NO 30 &lt;211&gt; LENGTH: 32 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 30	
taggatggtg atatggttga tacaggcttt gg	32
<pre>&lt;210&gt; SEQ ID NO 31 &lt;211&gt; LENGTH: 32 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (27)(27) &lt;223&gt; OTHER INFORMATION: I &lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (29)(29) &lt;223&gt; OTHER INFORMATION: I</pre>	
<400> SEQUENCE: 31	
taggatggtg atatggttga tacaggcttt gg	32
<pre>&lt;210&gt; SEQ ID NO 32 &lt;211&gt; LENGTH: 32 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (24)(24) &lt;223&gt; OTHER INFORMATION: I &lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (27)(27) &lt;223&gt; OTHER INFORMATION: I &lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (29)(29) &lt;223&gt; OTHER INFORMATION: I</pre>	
<400> SEQUENCE: 32	
taggatggtg atatggttga tacaggcttt gg	32
<210> SEQ ID NO 33 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 33	
tagggagact ttttgctaaa atgac	25

-co				

<210> SEQ ID NO 34 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 34 tagggatttt tgacccatct ttttctca 28 <210> SEQ ID NO 35 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 35 taggggccca accccatttt cat 23 <210> SEQ ID NO 36 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 36 25 taggtagaat gtttgcaatg caacc <210> SEQ ID NO 37 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 37 taggtcgcgt aatttgggta aggtcatc 28 <210> SEQ ID NO 38 <211> LENGTH: 30 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 38 tagtctcaaa gagaaagaga tcaaacaaga 30 <210> SEQ ID NO 39 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 39 tagtgcaatc catctagcag ctgt 24 <210> SEQ ID NO 40 <211> LENGTH: 29 <212> TYPE: DNA

-continued
------------

	-continued
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 40	
tataatccac ctatcccagt aggagaaat	29
<210> SEQ ID NO 41 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 41	
tatagggtca tgaatcaaga acccgg	26
<pre>&lt;210&gt; SEQ ID NO 42 &lt;211&gt; LENGTH: 26 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (24)(24) &lt;223&gt; OTHER INFORMATION: I</pre>	
<400> SEQUENCE: 42	
tatagggtca tgaatcaaga acccgg	26
<210> SEQ ID NO 43 <211> LENGTH: 32 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 43	
tatatatcac gagatctgca gtttatgagt aa	32
<210> SEQ ID NO 44 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 44	
tatatcacgt gatctgcagt ttatgag	27
<210> SEQ ID NO 45 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 45	
tatcagtgca atccatctag cagctgt	27
<210> SEQ ID NO 46 <211> LENGTH: 21 <212> TYPE: DNA	

<213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 46 21 tatcgctgga tgtgtctgcg g <210> SEQ ID NO 47 <211> LENGTH: 29 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 47 tatcttatca tgtctgggtc cccaggaag 29 <210> SEQ ID NO 48 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (18)...(18) <223> OTHER INFORMATION: I <400> SEQUENCE: 48 26 tatggtgcag tgggcatttg ataatg <210> SEQ ID NO 49 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 49 tatggtgcag tgggcatttg ataatg 26 <210> SEQ ID NO 50 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 50 tatgtcaaaa gctgtggaca atgat 25 <210> SEQ ID NO 51 <211> LENGTH: 18

18

<212> TYPE: DNA
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Primer

<400> SEQUENCE: 51

tatgtttgct atgcaacc

<210> SEQ ID NO 52 <211> LENGTH: 25 <212> TYPE: DNA 48

	-continued
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 52	
tcaaccacat aagaggtggg tgttc	25
<210> SEQ ID NO 53	
<211> LENGTH: 23 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 53	
tcaacctcca atcactcacc aac	23
<210> SEQ ID NO 54	
<211> LENGTH: 20 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 54	
tcaagagatc ttcagtttat	20
<210> SEQ ID NO 55	
<211> LENGTH: 26 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 55	
tcaagagatc ttcagtttat gagtaa	26
<210> SEQ ID NO 56	
<211> LENGTH: 26	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 56	
tcaagatggg gaatgagatt gccctt	26
cougaryyy guargagare geoere	20
<210> SEQ ID NO 57	
<211> LENGTH: 29 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 57	
tcaatgtatc ttatcatgtc tgggtcccc	29
<210> SEQ ID NO 58	
<211> LENGTH: 23	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	

-continued
------------

	-continued
<400> SEQUENCE: 58	
tcaccaaggt gcctctaaaa cga	23
<210> SEQ ID NO 59	
<211> LENGTH: 25	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 59	
tcaccacgca atgcttccct atctt	25
<210> SEQ ID NO 60	
<211> LENGTH: 24	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 60	
tcacctggga ccccatcgat ggac	24
<210> SEQ ID NO 61	
<211> LENGTH: 20	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 61	
tcacctttgc aaaggggccc	20
<210> SEQ ID NO 62	
<211> LENGTH: 26	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 62	
tcacgagatc tgcagtttat gagtaa	26
<210> SEQ ID NO 63	
<211> LENGTH: 25	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 63	
tcagagcatc agatcacctg ggacc	25
<210> SEQ ID NO 64	
<211> LENGTH: 29	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 64	
	29
tcaggatggt ttttggtaga ggctatagt	23

-continued

51

<210> SEQ ID NO 65 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 65 tcaggcaaga gtcactgcta tcgaga 26 <210> SEQ ID NO 66 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 66 tcagtgaata tgactagcca ggtgct 26 <210> SEQ ID NO 67 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 67 27 tcagtgtagg tagaatgttt gcaatgc <210> SEQ ID NO 68 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 68 tcagttgatg ggcctacctc atttctt 27 <210> SEQ ID NO 69 <211> LENGTH: 30 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 69 tcatctacaa catgatggga aagagagaga 30 <210> SEQ ID NO 70 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 70 tcatgtctgg gtcccctgga ag 22 <210> SEQ ID NO 71 <211> LENGTH: 29 <212> TYPE: DNA

	-continued
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 71	
~ tccaaccaga aatttgagaa cattgcaga	29
<210> SEQ ID NO 72 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 72	
tccacctgtg gttattgaac ctgt	24
<pre>&lt;210&gt; SEQ ID NO 73 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 73	
tccaggtact tcacccccaa acctg	25
<210> SEQ ID NO 74 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 74	
tccatacttc gccttcacca ttccc	25
<210> SEQ ID NO 75 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 75	
teccacatet gcaateggag acagg	25
<pre>&lt;210&gt; SEQ ID NO 76 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 76	
tccccaatct ccaatcactc accaa	25
<pre>&lt;210&gt; SEQ ID NO 77 &lt;211&gt; LENGTH: 24 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	

<400> SEQUENCE: 77	
tccccatcag acacctcagg tgga	24
<210> SEQ ID NO 78 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 78	
teeeeggtga ttggtateet ee	22
<210> SEQ ID NO 79 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 79	
tccctctaca gtagcaacgg atgcaa	26
<210> SEQ ID NO 80 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 80	
teegatgeae atteaeggtt ggtat	25
<210> SEQ ID NO 81 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 81	
tccgccttcg atgacagacg g	21
<pre>&lt;210&gt; SEQ ID NO 82 &lt;211&gt; LENGTH: 18 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 82	
tccggccctg gatgacaa	18
<210> SEQ ID NO 83 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 83	
tccgggctgt caatatgcta aaacg	25

-co				

<210> SEQ ID NO 84 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 84 tccgtcgcaa gttggtgatg agta 24 <210> SEQ ID NO 85 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 85 tcctaagagt gggaccatgg ataaa 25 <210> SEQ ID NO 86 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 86 28 tcctaagagt gggaccatgg ataaacac <210> SEQ ID NO 87 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 87 tcctagaagt agaggcagca tcca 24 <210> SEQ ID NO 88 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 88 teeteecaca geaaacatgt g 21 <210> SEQ ID NO 89 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 89 tcctcgggaa attggcattg cgat 24 <210> SEQ ID NO 90 <211> LENGTH: 24 <212> TYPE: DNA

	-continued
213> ORGANISM: Artificial Sequence	
220> FEATURE:	
223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 90	
contoteact cateagagea ecce	24
<210> SEQ ID NO 91	
<211> LENGTH: 22	
212> TYPE: DNA	
213> ORGANISM: Artificial Sequence	
220> FEATURE: 223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 91	
cctgcccac ctacaacaac ca	22
210> SEQ ID NO 92	
<211> LENGTH: 26	
(212> TYPE: DNA	
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	
223> OTHER INFORMATION: Primer	
<pre>&lt;400&gt; SEQUENCE: 92</pre>	
400> SEQUENCE: 92	
ccttcacgg aggctatgac taggta	26
<210> SEQ ID NO 93	
<211> LENGTH: 26	
<212> TYPE: DNA	
213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 93	
ccgaaaaacc tccaggcaag agtcac	26
<210> SEQ ID NO 94	
<211> LENGTH: 22	
<212> TYPE: DNA	
213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 94	
ccgcaagact gctagccgag ta	22
210. SEO TO NO 95	
<210> SEQ ID NO 95 <211> LENGTH: 25	
<212> TYPE: DNA	
213> ORGANISM: Artificial Sequence	
220> FEATURE:	
223> OTHER INFORMATION: Primer	
400> SEQUENCE: 95	
cgcagccat cttggcaaca tcctc	25
<210> SEQ ID NO 96	
211> LENGTH: 25 212> TYPE: DNA	
<211> LENGTH: 25	
2211> LENGTH: 25 2212> TYPE: DNA	

<400> SEQUENCE: 96		
togcagtooc caatotocaa toact	25	
<pre>&lt;210&gt; SEQ ID NO 97 &lt;211&gt; LENGTH: 20 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 97		
tcggcgcaag tcccagaacc	20	
<210> SEQ ID NO 98 <211> LENGTH: 20 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 98		
tcgggatgta atggctggtt	20	
<pre>&lt;210&gt; SEQ ID NO 99 &lt;211&gt; LENGTH: 26 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>		
<400> SEQUENCE: 99		
- tcgttatcgg ctcagctcta cagttc	26	
<pre>&lt;210&gt; SEQ ID NO 100 &lt;211&gt; LENGTH: 22 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>		
<400> SEQUENCE: 100		
tctacagtag caacggatgc aa	22	
<210> SEQ ID NO 101 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 101		
tctagccatg gcgttagtat gagtgt	26	
<210> SEQ ID NO 102 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 102		
tctcccgcac cctcaacaag c	21	

ued

<210> SEQ ID NO 103 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 103 tettttecag acceeggaet ee 22 <210> SEQ ID NO 104 <211> LENGTH: 19 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 104 tgaacatacc aatgcagtt 19 <210> SEQ ID NO 105 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 105 25 tgaactaacg gacagtggat atggc <210> SEQ ID NO 106 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 106 tgaatggttc cgccagtgca c 21 <210> SEQ ID NO 107 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 107 tgacacaggc acccacggaa ta 22 <210> SEQ ID NO 108 <211> LENGTH: 20 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 108 tgacgaacca cagcgtcaca 20 <210> SEQ ID NO 109 <211> LENGTH: 21 <212> TYPE: DNA

	-continued
<pre></pre> <213> ORGANISM: Artificial Sequence<220> FEATURE:<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 109	
tgacgtagcc aaggaggcgt t	21
<210> SEQ ID NO 110 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 110	
tgagaccatc ataagtagca agatgt	26
<210> SEQ ID NO 111 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 111	
tgagcagcag gaagcactat gg	22
<pre>&lt;210&gt; SEQ ID NO 112 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 112	
tgaggcctat agcagctata gcctc	25
<pre>&lt;210&gt; SEQ ID NO 113 &lt;211&gt; LENGTH: 24 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 113	
tgatgaaata tcgtggatgg aagc	24
<pre>&lt;210&gt; SEQ ID NO 114 &lt;211&gt; LENGTH: 24 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 114	
tgatgagata tcgtggatgg aagc	24
<210> SEQ ID NO 115 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	

<400> SEQUENCE: 115	
tgattgaccc ttttcagttg ggcctt	26
<pre>&lt;210&gt; SEQ ID NO 116 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 116	
tgattgatgt tgaaaccacc gaaga	25
<pre>&lt;210&gt; SEQ ID NO 117 &lt;211&gt; LENGTH: 24 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 117</pre>	
<400> SEQUENCE: 117	24
tgcagtcccc aatctccaat cact <210> SEQ ID NO 118 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	24
<400> SEQUENCE: 118	
tgccaatcac tcatacaacc cccaa	25
<pre>&lt;210&gt; SEQ ID NO 119 &lt;211&gt; LENGTH: 24 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 119	
tgccacgacc ggcaagacca acat	24
<pre>&lt;210&gt; SEQ ID NO 120 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 120	
tgccagatga tcctataatt gttgagga	28
<pre>&lt;210&gt; SEQ ID NO 121 &lt;211&gt; LENGTH: 26 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 121	
tgccagggag tagtagaagc aatgaa	26

-co				

<210> SEQ ID NO 122 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 122 tgccattcat aggctgccca tc 22 <210> SEQ ID NO 123 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 123 tgcctacatt ctatcttgcg ttacatgg 28 <210> SEQ ID NO 124 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 124 28 tgcctagtgt gaagatgggg aatgagat <210> SEQ ID NO 125 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 125 tgcctgaatt ggagatatga gaccat 26 <210> SEQ ID NO 126 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 126 tgccttaatt ggagatatga gac 23 <210> SEQ ID NO 127 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 127 tgccttaatt ggagatatga gaccat 26 <210> SEQ ID NO 128 <211> LENGTH: 29 <212> TYPE: DNA

	-continued
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 128	
tgcgttaact ggaccaatga gaactttcc	29
<210> SEQ ID NO 129 <211> LENGTH: 25	
<pre>&lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence</pre>	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 129	25
tgctcctctc aaaaccatgg gacac	25
<210> SEQ ID NO 130 <211> LENGTH: 23	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 130	
tgctgaattc gtaacattga tga	23
<210> SEQ ID NO 131 <211> LENGTH: 28	
<212> TYPE: DNA	
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 131	
tggaaateet tttteteaag gaegtggt	28
<210> SEQ ID NO 132	
<211> LENGTH: 23	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 132	
tggacattca tctctatcag tgc	23
<210> SEQ ID NO 133	
<211> LENGTH: 22 <212> TYPE: DNA	
<pre><li><li><li><li><li><li><li><li><li><li< td=""><td></td></li<></li></li></li></li></li></li></li></li></li></pre>	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 133	
tggaccatca ttggaaggga cg	22
<210> SEQ ID NO 134	
<211> LENGTH: 24 <212> TYPE: DNA	
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	
<223> OTHER INFORMATION: Primer	

<400> SEQUENCE: 134		
tggaggctcc gttgtctgca tgta	24	
<210> SEQ ID NO 135		
<211> LENGTH: 30		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 135		
tggatagagg agaatgaatg gatggaagac	30	
<210> SEQ ID NO 136		
<211> LENGTH: 19		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 136		
tggatgtgtc tgcggcgtt	19	
<210> SEQ ID NO 137		
<211> SEQ 1D NO 137 <211> LENGTH: 19		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 137		
tggccccgtg tctggtagc	19	
	19	
<210> SEQ ID NO 138	19	
<210> SEQ ID NO 138 <211> LENGTH: 28	19	
<210> SEQ ID NO 138 <211> LENGTH: 28 <212> TYPE: DNA	19	
<210> SEQ ID NO 138 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence	19	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	19	
<210> SEQ ID NO 138 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence	19	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	19	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	19 28	
<210> SEQ ID NO 138 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 138		
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139</pre>		
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25</pre>		
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA</pre>		
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence</pre>		
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA</pre>		
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>		
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139</pre>	28	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>		
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139 tggctgaact tttgaaagtg agcgg &lt;210&gt; SEQ ID NO 140</pre>	28	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139 tggctgaact tttgaaagtg agcgg &lt;210&gt; SEQ ID NO 140 &lt;211&gt; LENGTH: 25</pre>	28	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139 tggctgaact tttgaaagtg agcgg &lt;210&gt; SEQ ID NO 140 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA</pre>	28	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139 tggctgaact tttgaaagtg agcgg &lt;210&gt; SEQ ID NO 140 &lt;211&gt; LENGTH: 25</pre>	28	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139 tggctgaact tttgaaagtg agcgg &lt;210&gt; SEQ ID NO 140 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA</pre>	28	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139 tggctgaact tttgaaagtg agcgg &lt;210&gt; SEQ ID NO 140 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence</pre>	28	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139 tggctgaact tttgaaagtg agcgg &lt;210&gt; SEQ ID NO 140 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;20&gt; FEATURE: &lt;220&gt; SEQ ID NO 140 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	28	
<pre>&lt;210&gt; SEQ ID NO 138 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 138 tggctaggtg tcacgtagaa aggactac &lt;210&gt; SEQ ID NO 139 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 139 tggctgaact tttgaaagtg agcgg &lt;210&gt; SEQ ID NO 140 &lt;211&gt; LENGTH: 25 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	28	

-co			

<210> SEQ ID NO 141 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 141 tgggatgatc tactgagact gtgtga 26 <210> SEQ ID NO 142 <211> LENGTH: 30 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 142 tgggggtatt ttgacttcaa ccgattccac 30 <210> SEQ ID NO 143 <211> LENGTH: 30 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (22)...(22) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (24)...(24) <223> OTHER INFORMATION: I <400> SEQUENCE: 143 30 tgggggtatt ttgacttcaa ccgattccac <210> SEQ ID NO 144 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 144 tgggtcctgc ccacctacaa caacca 26 <210> SEQ ID NO 145 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (4)...(4) <223> OTHER INFORMATION: I <400> SEQUENCE: 145 tggggtggtg agagtttggg a 21 <210> SEQ ID NO 146 <211> LENGTH: 22 <212> TYPE: DNA

	-continued
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Drimor</pre>	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 146	20
tggtacteet caactgetgt ee	22
<210> SEQ ID NO 147 <211> LENGTH: 18 <212> TYPE: DNA	
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 147	
tggtgccgac ggagtggg	18
<210> SEQ ID NO 148 <211> LENGTH: 18 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 148	
tggtgccgat ggagtggg	18
<pre>&lt;210&gt; SEQ ID NO 149 &lt;211&gt; LENGTH: 21 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 149</pre>	
tggtgtggtg agagtttggg a	21
<210> SEQ ID NO 150 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 150	
tggtgttgtg agggtctggg a	21
<210> SEQ ID NO 151 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 151	
tgtaaatcta gtggctctcc tcccac	26
<pre>&lt;210&gt; SEQ ID NO 152 &lt;211&gt; LENGTH: 32 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	

<400> SEQUENCE: 152	
tgtacaatct atgtaggaga tccttactgt cc	32
210, CEO TO NO 152	
<210> SEQ ID NO 153	
<211> LENGTH: 23	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 153	
tgtaggtaga atgtttgcaa tgc	23
<210> SEQ ID NO 154	
<211> LENGTH: 24	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 154	
tgtaggtgat cccttcaatc ctcc	24
<210> SEQ ID NO 155	
<211> LENGTH: 28	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 155	
tgtcacactt gcatttacaa catgatgg	28
<210> SEQ ID NO 156	
<211> LENGTH: 21	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
(223) OTHER INFORMATION. FILMEI	
<400> SEQUENCE: 156	
tgtcccgact gtggcacttg a	21
<210> SEQ ID NO 157	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 157	
tgtctatctt cctgactc	18
<210> SEQ ID NO 158	
<211> LENGTH: 23	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<220> FEATURE:	
<221> NAME/KEY: modified_base	
<222> LOCATION: (9)(9)	
<223> OTHER INFORMATION: I	
<220> FEATURE:	

-cor	nti	.nu	ιed

<221> NAME/KEY: modified\_base <222> LOCATION: (12)...(12) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (18)...(18) <223> OTHER INFORMATION: I <400> SEQUENCE: 158 tgtctctctg gcgatgggtg gcc 23 <210> SEQ ID NO 159 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 159 tgtctctctg gagatgggtg gcc 23 <210> SEQ ID NO 160 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 160 23 tgtctttcag gagatggatg gcc <210> SEQ ID NO 161 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (12)...(12) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (15)...(15) <223> OTHER INFORMATION: I <400> SEQUENCE: 161 tgtctttcag gagatggatg gcc 23 <210> SEQ ID NO 162 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 162 tgtgaagcgc acaattgtta agac 24 <210> SEQ ID NO 163 <211> LENGTH: 32 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer

<400> SEQUENCE: 163	
tgtgaatatc aagcaggaca taacaaggta gg	32
<210> SEQ ID NO 164 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 164	
tgtgcatctt actcactaaa ggagaa	26
eyeyeaceee accounced yyayaa	23
<pre>&lt;210&gt; SEQ ID NO 165 &lt;211&gt; LENGTH: 24 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 165	
tgtgccctca aaaaccctaa cctc	24
<210> SEQ ID NO 166 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 166	
tgtgccctca aaaaccctga cctc	24
<210> SEQ ID NO 167 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 167	
tgtggacaat gatctccatt gctgcaat	28
<210> SEQ ID NO 168 <211> LENGTH: 20 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 168	
tgtgtcatct gcgagtgtgg	20
<210> SEQ ID NO 169 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 169	
tgttaatgtc tatcttcctg actc	24

-continued

<210> SEQ ID NO 170 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 170 tgtttgcaca gaggctaaat ga 22 <210> SEQ ID NO 171 <211> LENGTH: 29 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 171 ttacaaagag ggccaacgcc attctcatc 29 <210> SEQ ID NO 172 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 172 25 ttatgtcctt ggggtcacct actgg <210> SEQ ID NO 173 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 173 ttcaatcaag catttcggta tgaaac 26 <210> SEQ ID NO 174 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 174 ttcacgcaga aagcgtctag cca 23 <210> SEQ ID NO 175 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 175 ttcacgcaga aagcgtctag ccat 24 <210> SEQ ID NO 176 <211> LENGTH: 24 <212> TYPE: DNA

	-continued
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 176	
ttcagatgtc tgtgtggcgg ccta	24
<210> SEQ ID NO 177	
<211> LENGTH: 24	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<220> FEATURE: <221> NAME/KEY: modified_base	
<2221> NAME/RET: Modified_base <222> LOCATION: (19)(19)	
<223> OTHER INFORMATION: I	
<400> SEQUENCE: 177	
ttcagatgtc tgtgtggcgg ccta	24
<210> SEQ ID NO 178 <211> LENGTH: 24	
<211> LENGTH: 24 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<220> FEATURE:	
<221> NAME/KEY: modified_base	
<222> LOCATION: (12)(12) <223> OTHER INFORMATION: I	
<220> FEATURE:	
<221> NAME/KEY: modified_base	
<222> LOCATION: (19)(19) <223> OTHER INFORMATION: I	
<400> SEQUENCE: 178	
ttcagatgtc tgtgtggcgg ccta	24
<210> SEQ ID NO 179	
<211> LENGTH: 28	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 179	
ttgcaaagag gcccaacccc attctcat	28
<210> SEQ ID NO 180	
<211> LENGTH: 28	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<220> FEATURE: <221> NAME/KEY: modified_base	
<pre>&lt;222&gt; LOCATION: (24)(24) &lt;223&gt; OTHER INFORMATION: I</pre>	
<400> SEQUENCE: 180	
ttgcaaagag gcccaacccc attctcat	28
<210> SEQ ID NO 181	
<211> LENGTH: 25 <212> TYPE: DNA	
-continued	
------------	
------------	

	-continued
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 181	ar.
ttgcgtcaat ggcttatatc aaagg	25
<210> SEQ ID NO 182 <211> LENGTH: 25	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 182	25
ttggatcagc cactgatgaa agatc	25
<210> SEQ ID NO 183	
<211> LENGTH: 24 <212> TYPE: DNA	
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	
<2205 FLATORE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 183	
ttggtccatg actccgacct tgaa	24
<210> SEQ ID NO 184	
<211> LENGTH: 24	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 184	
tttcgttgtg gcaaggtagg acac	24
<210> SEQ ID NO 185	
<211> LENGTH: 17	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 185	17
tttggacacc caatggt	17
<210> SEQ ID NO 186	
<211> LENGTH: 23 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 186	
taaggaacgg cgacaggatg cgc	23
<210> SEQ ID NO 187	
<211> LENGTH: 30	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	

=0	continued
<400> SEQUENCE: 187	
taatagactt tcccctctat attctaattc	30
-	
<210> SEQ ID NO 188	
<211> LENGTH: 30	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 188	
taataggett teceetetat attetaatte	30
010. CTO TO NO 100	
<210> SEQ ID NO 189 <211> LENGTH: 30	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 189	
taatccagaa gaaccaacaa gaagatgagg	30
<210> SEQ ID NO 190	
<211> LENGTH: 19	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 190	
taattgcatt agtctcact	19
<210> SEQ ID NO 191	
<211> LENGTH: 16	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 191	
tacagaatag gctttg	16
<210> SEQ ID NO 192	
<211> LENGTH: 30 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 192	
tacatattca tccgtgctta caaccttaga	30
<210> SEQ ID NO 193	
<211> LENGTH: 24	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<2220> FLATORE: <2223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 193	
taccatgtcc gaacctgtat ctgt	24

	continued
<210> SEQ ID NO 194	
<211> LENGTH: 21	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 194	
tacggagagt gacgggtttc c	21
<210> SEQ ID NO 195	
<211> LENGTH: 25	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 195	
tacgggacgt aaacaaagga cgtcc	25
<210> SEQ ID NO 196	
<211> LENGTH: 20	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 196	
tacggggcca tgccgtgttg	20
<210> SEQ ID NO 197	
<211> LENGTH: 26	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 197	
tactcaccgg ttccgcagac cactat	26
<210> SEQ ID NO 198	
<211> LENGTH: 26	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 198	
tactccgtct cgtacgactt tctgtt	26
<210> SEQ ID NO 199	
<211> LENGTH: 18	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 199	
tactctaaca gcatccat	18
2	
<210> SEQ ID NO 200	
<211> LENGTH: 30	
<212> TYPE: DNA	

-continued		
<213> ORGANISM: Artificial Sequence <220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 200		
tactcttgtt gtcacagcta tagcttgatt	30	
<210> SEQ ID NO 201 <211> LENGTH: 34 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 201		
tactgcataa tagactttcc cctctatatt ctaa	34	
<210> SEQ ID NO 202 <211> LENGTH: 34 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 202		
tactgcataa taggctttcc cctctaaatt ctaa	34	
<210> SEQ ID NO 203 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 203		
tacttgggat tgtttgtgtg agacgg	26	
<210> SEQ ID NO 204 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 204		
tagatttcat ttagcctctg tgcaaa	26	
<210> SEQ ID NO 205 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 205		
tagcaaattg ctttgctgtt gcacta	26	
<210> SEQ ID NO 206 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer		

	-continued
<400> SEQUENCE: 206	
tageteegag ceacatgaae e	21
<210> SEQ ID NO 207	
<211> LENGTH: 24 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 207	
tagctgaaat tgctgctgga gagg	24
<210> SEQ ID NO 208	
<211> LENGTH: 28	
<pre>&lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence</pre>	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 208	
taggaatatg ataaaacgcc gcagacac	28
<210> SEQ ID NO 209	
<211> LENGTH: 25	
<212> TYPE: DNA	
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 209	
tagggagatt gggageteee gtatt	25
<210> SEQ ID NO 210	
<211> LENGTH: 25	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<pre>&lt;220&gt; FEATURE:</pre>	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 210	
tagggtcact tggaggattg aaagg	25
<210> SEQ ID NO 211	
<211> LENGTH: 30	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer <400> SEQUENCE: 211	
tagtattttg tcctgccacg catttaaacg	30
<210> SEQ ID NO 212	
<211> LENGTH: 30 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<pre>&lt;223&gt; OTHER INFORMATION: Primer &lt;220&gt; FEATURE:</pre>	
<221> NAME/KEY: modified_base	
<222> LOCATION: (20)(20)	
<223> OTHER INFORMATION: I	

-continued		
<400> SEQUENCE: 212		
tagtattttg tcctgccacg catttaaacg	30	
<210> SEQ ID NO 213		
<211> LENGTH: 30 <212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer <220> FEATURE:		
<221> NAME/KEY: modified_base		
<222> LOCATION: (19)(20) <223> OTHER INFORMATION: I		
<400> SEQUENCE: 213		
tagtattttg tcctgccacg catttaaacg	30	
<210> SEQ ID NO 214		
<211> LENGTH: 30 <212> TYPE: DNA		
<212> TIPE: DNA <213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer <220> FEATURE:		
<220> FLAIORE: <221> NAME/KEY: modified_base		
<222> LOCATION: (18)(20)		
<223> OTHER INFORMATION: I		
<400> SEQUENCE: 214		
tagtattttg teetgecaeg catttaaaeg	30	
<210> SEQ ID NO 215		
<211> LENGTH: 25		
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 215		
tagtcgttgg tccgttgtta gggaa	25	
<210> SEQ ID NO 216		
<211> LENGTH: 24 <212> TYPE: DNA		
<212> TIPE: DNA <213> ORGANISM: Artificial Sequence		
<220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 216		
- tatatgataa aacgccgcag acac	24	
<210> SEQ ID NO 217		
<211> LENGTH: 25 <212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 217		
tatccccgaa agcccagata tatac	25	
-210, SEO TO NO 219		
<210> SEQ ID NO 218 <211> LENGTH: 28		
<212> TYPE: DNA		

-	-continued
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 218	
tatttgcctg ccaattgctt tttaaaaa	28
<210> SEQ ID NO 219	
<211> LENGTH: 28 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<220> FEATURE:	
<221> NAME/KEY: modified_base	
<222> LOCATION: (12)(13) <223> OTHER INFORMATION: I	
<220> FEATURE:	
<221> NAME/KEY: modified_base <222> LOCATION: (20)(20)	
<223> OTHER INFORMATION: I	
<400> SEQUENCE: 219	
tatttgcctg ccaattgctt tttaaaaa	28
<210> SEQ ID NO 220	
<211> LENGTH: 23	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 220	
tcaaacccag aggtgcctgt aaa	23
<210> SEQ ID NO 221	
<211> LENGTH: 20 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 221	
tcaacactgg tgttgtccca	20
	20
<210> SEQ ID NO 222	
<211> LENGTH: 26 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 222	
tcaatctctc ccttaaccat cccatt	26
<210> SEQ ID NO 223	
<211> LENGTH: 25	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 223	
tcactcgttt agacggcgga atatc	25

<210> SEQ ID NO 224 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 224 tcagcaaatt gttctgctgc tgcacta 27 <210> SEQ ID NO 225 <211> LENGTH: 29 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 225 tcagctatca atttctctgc caatatttg 29 <210> SEQ ID NO 226 <211> LENGTH: 29 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 226 29 tcagctatca ttttctcagc caagatttg <210> SEQ ID NO 227 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 227 tcaggagctt tcgtgtccac ctt 23 <210> SEQ ID NO 228 <211> LENGTH: 19 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 228 tcaggcacag cttggaggc 19 <210> SEQ ID NO 229 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 229 tcaggctgcc acaccagatg tc 22 <210> SEQ ID NO 230 <211> LENGTH: 17 <212> TYPE: DNA

-continu	ed
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 230	
tcaggggtcc ttctata	17
<210> SEQ ID NO 231	
<211> LENGTH: 20	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 231	
ccaggggtgc tttgtgcttc	20
<210> SEQ ID NO 232 <211> LENGTH: 22	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 232	
tcaggtggtg tctgccctcg tt	22
<210> SEQ ID NO 233	
<211> LENGTH: 22 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 233	
tcattgtggt gggtaggtcg tc	22
<210> SEQ ID NO 234 <211> LENGTH: 20	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence <220> FEATURE:	
<pre>&lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 234	
tcatttagcc tctgtgcaaa	20
<210> SEQ ID NO 235	
<211> LENGTH: 25	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 235	
tccacatgaa ccaaatggct ctgct	25
<210> SEQ ID NO 236	
<211> LENGTH: 28	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	

-continued			
<400> SEQUENCE: 236			
tccagaagaa ccaacaagaa gatgaggc	28		
<210> SEQ ID NO 237 <211> LENGTH: 33			
<212> TYPE: DNA			
<213> ORGANISM: Artificial Sequence <220> FEATURE:			
<223> OTHER INFORMATION: Primer			
<pre>&lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: modified_base</pre>			
<222> LOCATION: (11)(11)			
<223> OTHER INFORMATION: I			
<400> SEQUENCE: 237			
tccagaccca cccaagaatg agaacacaag ata	33		
<210> SEQ ID NO 238			
<211> LENGTH: 27 <212> TYPE: DNA			
<213> ORGANISM: Artificial Sequence			
<220> FEATURE: <223> OTHER INFORMATION: Primer			
<400> SEQUENCE: 238			
tccagaccct gcaaaaaatg agaacac	27		
<210> SEQ ID NO 239			
<211> LENGTH: 24			
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence			
<220> FEATURE: <223> OTHER INFORMATION: Primer			
<400> SEQUENCE: 239			
tccatatacg aacacccgg cgac	24		
<210> SEQ ID NO 240			
<211> LENGTH: 21 <212> TYPE: DNA			
<213> ORGANISM: Artificial Sequence			
<220> FEATURE: <223> OTHER INFORMATION: Primer			
<400> SEQUENCE: 240			
tccatcgatg gggtcccagg t	21		
<210> SEQ ID NO 241 <211> LENGTH: 24			
<212> TYPE: DNA			
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>			
<223> OTHER INFORMATION: Primer			
<400> SEQUENCE: 241			
tccatgggga atcgcaatgc caat	24		
<210> SEQ ID NO 242			
<211> LENGTH: 24 <212> TYPE: DNA			
<213> ORGANISM: Artificial Sequence			
<220> FEATURE: <223> OTHER INFORMATION: Primer			

			u	

	-continued
<400> SEQUENCE: 242	
tccatttgga atcgcaatgc caat	24
<210> SEQ ID NO 243	
<211> LENGTH: 22	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 243	
tccgatgtgg cttatttggc tg	22
<210> SEQ ID NO 244	
<211> LENGTH: 27	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 244	
tccgtggatt gttctgtagc agtcttc	27
<210> SEQ ID NO 245	
<211> LENGTH: 32	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 245	
tectatteet ceettettt taaaatteat ge	32
<210> SEQ ID NO 246	
<211> LENGTH: 23	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 246	
teeteggtag teettetae gtg	23
<210> SEQ ID NO 247	
<211> LENGTH: 24	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 247	
teetgettga tggeetgtaa gtet	24
<210> SEQ ID NO 248	
<211> LENGTH: 20	
<212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 248	
tcctttgcca tcccattgtc	20

	-continued
<210> SEQ ID NO 249 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 249	
tcgacacatt ggaggagcat gatgtta	27
<210> SEQ ID NO 250 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 250	
tcgatcaggc tcaacatagc cctctt	26
<pre>&lt;210&gt; SEQ ID NO 251 &lt;211&gt; LENGTH: 22 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 251</pre>	
tegcaageac ectateagge ag	22
<210> SEQ ID NO 252 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 252	
tcgccatgtc tctctacctg cg	22
<210> SEQ ID NO 253 <211> LENGTH: 32 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 253	
tcggcagcca tttgcaaata atcaggatat tt	32
<210> SEQ ID NO 254 <211> LENGTH: 20 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 254	
tcgggagggc aatgagggtc	20
<210> SEQ ID NO 255 <211> LENGTH: 23 <212> TYPE: DNA	

	-continued
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 255	
	23
tcgggcccac tatgacgtgt aca	23
<210> SEQ ID NO 256	
<211> LENGTH: 20 <212> TYPE: DNA	
<212> IIPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 256	
tcgggtactg gtgcacacgc	20
<210> SEQ ID NO 257	
<211> LENGTH: 28	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 257	
tcgtttgtag ggaacattgg tgaggaag	28
<210> SEQ ID NO 258	
<211> LENGTH: 28	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 258	
tcgtttgtag ggaacattgg tgaggaag	28
<210> SEQ ID NO 259	
<211> LENGTH: 25 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 259	
	25
tctaaacgca catccccatg aattt	25
<210> SEQ ID NO 260	
<211> LENGTH: 23	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 260	
tctaaatgag gacctgacct gtg	23
<210> SEQ ID NO 261	
<211> LENGTH: 20 <212> TYPE: DNA	
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence	
<220> FEATURE:	
<223> OTHER INFORMATION: Primer	

# -continued

<400> SEQUENCE: 261		
tctaatattg tgtttatcca	20	
-210- CEO TO NO 262		
<210> SEQ ID NO 262 <211> LENGTH: 26		
<211> LENGIR: 26 <212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 262		
tctatgagga cattttccca atctcc	26	
	20	
<210> SEQ ID NO 263		
<211> LENGTH: 27		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 263		
tctcccatca tagtatatcc ttttgct	27	
<210> SEQ ID NO 264		
<211> LENGTH: 24		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 264		
tctctgggtg gggaaggagg ggag	24	
<210> SEQ ID NO 265		
<211> LENGTH: 32		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 265		
tctgcaacca tttgaaaata atctggatat tt	32	
<210> SEQ ID NO 266		
<211> LENGTH: 32		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 266		
tctgcaacca tttgaaaata atctggatat tt	32	
<210> SEQ ID NO 267		
<211> LENGTH: 32		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 267		
tctgcaacca tttgcaaata atctggatat tt	32	

-continued

<210> SEO ID NO 268 <211> LENGTH: 35 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 268 tctgcaacca tttgcaaata atctggatat ttgca 35 <210> SEQ ID NO 269 <211> LENGTH: 35 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (33)...(33) <223> OTHER INFORMATION: I <400> SEQUENCE: 269 tctgcaacca tttgcaaata atctggatat ttgca 35 <210> SEQ ID NO 270 <211> LENGTH: 32 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (15)...(15) <223> OTHER INFORMATION: I <400> SEQUENCE: 270 32 tctgcaacca tttgcaaata atctggatat tt <210> SEQ ID NO 271 <211> LENGTH: 35 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (15)...(15) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (33)...(33) <223> OTHER INFORMATION: I <400> SEQUENCE: 271 tctgcaacca tttgcaaata atctggatat ttgca 35 <210> SEQ ID NO 272 <211> LENGTH: 32 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (14)...(15) <223> OTHER INFORMATION: I

-continued

<400> SEQUENCE: 272 32 tctgcaacca tttgcaaata atctggatat tt <210> SEQ ID NO 273 <211> LENGTH: 35 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (14)...(15) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (33)...(33) <223> OTHER INFORMATION: I <400> SEQUENCE: 273 tctgcaacca tttgcaaata atctggatat ttgca 35 <210> SEQ ID NO 274 <211> LENGTH: 32 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (14)...(15) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (17)...(17) <223> OTHER INFORMATION: I <400> SEQUENCE: 274 32 tctgcaacca tttgcaaata atctggatat tt <210> SEQ ID NO 275 <211> LENGTH: 19 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 275 tctggaaccc agcagtgga 19 <210> SEQ ID NO 276 <211> LENGTH: 29 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 276 tcttcatcag agctaggaga ttctgttgt 29 <210> SEQ ID NO 277 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer

# -continued

<400> SEQUENCE: 277		
(400) BEQUENCE. 277		
tcttgctttg acatgttgtg gtgga	25	
<210> SEQ ID NO 278		
<211> LENGTH: 21		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 278		
tctttcccct ctgtattcta a	21	
<210> SEQ ID NO 279		
<211> LENGTH: 37		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 279		
	37	
tgaaaacaca agatactttg gaacatacac aggaggt	57	
<210> SEQ ID NO 280		
<211> LENGTH: 24		
<212> TYPE: DNA		
<pre>&lt;213&gt; ORGANISM: Artificial Sequence</pre>		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 280		
tgaaccaatt gcattagtct cact	24	
<210> SEQ ID NO 281		
<211> LENGTH: 29		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 281		
tgaagaggaa tatgataaaa cgccgcaga	29	
-210- SEO ID NO 282		
<210> SEQ ID NO 282 <211> LENGTH: 28		
<211> LENGIN: 28 <212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE: <223> OTHER INFORMATION: Primer		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 282		
tgaagggaaa gttctcattg gtccagtt	28	
<210> SEQ ID NO 283		
<211> LENGTH: 28		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 283		
tgaagtcatc cagtatagtg tttatcca	28	
-jj-sate cajeacajej tetateta	20	

-	С	0	n	t	i	n	u	e	d

<210> SEQ ID NO 284 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 284 tgaatgtacc ccatgaggtc ggc 23 <210> SEQ ID NO 285 <211> LENGTH: 37 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (18)...(18) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (27)...(27) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (34)...(34) <223> OTHER INFORMATION: I <400> SEQUENCE: 285 37 tgagaacaca agatactttg gaaactttca caggtgg <210> SEQ ID NO 286 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 286 tgagagtgga tgggcagcct atg 23 <210> SEQ ID NO 287 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 287 tgaggatgtg tccaagatgg ctgcg 25 <210> SEQ ID NO 288 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 288 tgaggctgga tctatccacg caaa 24 <210> SEQ ID NO 289 <211> LENGTH: 25 <212> TYPE: DNA

-cont	tinued	
ORGANISM: Artificial Sequence FEATURE:		
OTHER INFORMATION: Primer		
SEQUENCE: 289		
ggagt cgagggataa ggaac	25	
SEQ ID NO 290		
LENGTH: 27		
TYPE: DNA		
ORGANISM: Artificial Sequence		
FEATURE: OTHER INFORMATION: Primer		
SEQUENCE: 290		
gggat aatggcatag gaggaat	27	
SEQ ID NO 291 LENGTH: 23		
TYPE: DNA		
ORGANISM: Artificial Sequence		
FEATURE:		
OTHER INFORMATION: Primer		
SEQUENCE: 291		
gattg gcggggtaag gac	23	
SEQ ID NO 292		
LENGTH: 23		
TYPE: DNA		
ORGANISM: Artificial Sequence FEATURE:		
OTHER INFORMATION: Primer		
SEQUENCE: 292		
aaacg ccgcagacac atc	23	
SEQ ID NO 293		
LENGTH: 23		
TYPE: DNA		
ORGANISM: Artificial Sequence FEATURE:		
OTHER INFORMATION: Primer		
SEQUENCE: 293		
gtcgc cttgaaccat tgc	23	
SEQ ID NO 294		
LENGTH: 21		
TYPE: DNA		
ORGANISM: Artificial Sequence FEATURE:		
OTHER INFORMATION: Primer		
SEQUENCE: 294		
acaac ggacacacaa a	21	
SEQ ID NO 295		
LENGTH: 33		
TYPE: DNA		
LENGTH: 33		

# -continued

<400> SEQUENCE: 295	
tgcactagaa aacttcatct gggttgcagt atc	33
<210> SEQ ID NO 296 <211> LENGTH: 33 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 296	
tgcactagag aatttcatct gggttgcagt atc	33
<pre>&lt;210&gt; SEQ ID NO 297 &lt;211&gt; LENGTH: 20 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer </pre>	
<400> SEQUENCE: 297	
<pre>c210&gt; SEQ ID NO 298 c211&gt; VENCEUL 22</pre>	20
<211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 298	
tgcatagcaa tctcaggagg tgc	23
<210> SEQ ID NO 299 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 299	
tgcatgaata agggtctgaa gccca	25
<210> SEQ ID NO 300 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 300	
tgcatgttga tgttgtccag catgat	26
<210> SEQ ID NO 301 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 301	
tgccgaggcc tatgtggtcg acatt	25

-co			

<210> SEQ ID NO 302 <211> LENGTH: 31 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 302 tgcctgtgct tccaaggaat tgtgtgtaat a 31 <210> SEQ ID NO 303 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 303 tgcgggaaag gcgggagttg aagat 25 <210> SEQ ID NO 304 <211> LENGTH: 18 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 304 18 tgctaacaat ttctcagc <210> SEQ ID NO 305 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 305 tgctccagga ggtgcaaatc aaaga 25 <210> SEQ ID NO 306 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 306 tgctccccct agaaaattga gagaagtc 28 <210> SEQ ID NO 307 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 307 tgctgatctt ggctgtccac ctc 23 <210> SEQ ID NO 308 <211> LENGTH: 23 <212> TYPE: DNA

	-continued
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 308	
- tgctggcatg aaccaaactg att	23
<210> SEQ ID NO 309 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 309	
tgcttccatc cacgatatct catc	24
<210> SEQ ID NO 310 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 310	
tggaccatac aggcaacccg acaa	24
<pre>&lt;210&gt; SEQ ID NO 311 &lt;211&gt; LENGTH: 22 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 311</pre>	
tggcaattcc ggtgtactca cc	22
<210> SEQ ID NO 312 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 312	
tggcagtgga atcggttgaa gtcaaa	26
<210> SEQ ID NO 313 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 313	
tggcatggtg ccaaactgat t	21
<210> SEQ ID NO 314 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE:	

-continued				
<pre>&lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (7)(7) &lt;223&gt; OTHER INFORMATION: I</pre>				
<400> SEQUENCE: 314				
tggcatggtg ccaaactgat t	21			
<210> SEQ ID NO 315 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer				
<400> SEQUENCE: 315				
tggcttgatt gtctccttga accattgc	28			
<210> SEQ ID NO 316 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer				
<400> SEQUENCE: 316				
tgggaaagct ggtagaggta catgc	25			
<210> SEQ ID NO 317 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer				
<400> SEQUENCE: 317				
tgggactctg gttcatgtat tgg	23			
<210> SEQ ID NO 318 <211> LENGTH: 18 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer				
<400> SEQUENCE: 318				
tgggcaatga gggtcact	18			
<210> SEQ ID NO 319 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer				
<400> SEQUENCE: 319				
tgggcctctc tgcaaaggag a	21			
<pre>&lt;210&gt; SEQ ID NO 320 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: </pre>				

<223> OTHER INFORMATION: Primer

# -continued

<400> SEQUENCE: 320	
tgggctgata tggaatgtgc cctatctg	28
<210> SEQ ID NO 321 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 321	
tggggagatg accatgatgg tgc	23
<pre>&lt;210&gt; SEQ ID NO 322 &lt;211&gt; LENGTH: 23 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 322</pre>	
tggggctcat ggtcattgtc atc	23
<210> SEQ ID NO 323 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 323	
tgggggacct aattgctact gtatctga	28
<210> SEQ ID NO 324 <211> LENGTH: 29 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 324	
tgggtccctg atccaactag aaatgaaaa	29
<210> SEQ ID NO 325 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 325	
tggtaatcaa aatactgcgg gccaa	25
<210> SEQ ID NO 326 <211> LENGTH: 33 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 326	
tggtactctt gatgtcacag ctatagcttg att	33

-continued

<210> SEQ ID NO 327 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 327 tggtcaaagc tggcatgaac caaa 24 <210> SEQ ID NO 328 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (7)...(7) <223> OTHER INFORMATION: I <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (11)...(11) <223> OTHER INFORMATION: I <400> SEQUENCE: 328 tggtcaaagc tggcatgaac caaa 24 <210> SEQ ID NO 329 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 329 24 tggtctggaa aagacagggt tggg <210> SEQ ID NO 330 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 330 24 tggtcttctg gggcttgttc catc <210> SEQ ID NO 331 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 331 tggtgtttgg agtggctatt ggcag 25 <210> SEQ ID NO 332 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer

# -continued

<400> SEQUENCE: 332	
tgtaaccagc agaataggct ttg	23
<210> SEQ ID NO 333 <211> LENGTH: 31 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 333	
tgtaactgtt aaaactgagg ttgttggagt g	31
<pre>&lt;210&gt; SEQ ID NO 334 &lt;211&gt; LENGTH: 22 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 334</pre>	
tgtacggcag tgtttgcggt at	22
<210> SEQ ID NO 335 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 335	25
tgtagaggtg gttgttgtag gtggg	25
<210> SEQ ID NO 336 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 336	
tgtatctgaa gctgctgctg c	21
<pre>&lt;210&gt; SEQ ID NO 337 &lt;211&gt; LENGTH: 28 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 337	
tgtcatatcc cctattcctc cccttctt	28
<pre>&lt;210&gt; SEQ ID NO 338 &lt;211&gt; LENGTH: 24 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 338	
tgtcatcttg ccctcctccc acca	24

-c			

<210> SEQ ID NO 339 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 339 tgtcattttg ccctcctccc acca 24 <210> SEQ ID NO 340 <211> LENGTH: 30 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 340 tgtccataat ttctggcaat aaccttctat 30 <210> SEQ ID NO 341 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 341 27 tgtctaaacg cacatcccca tgaattt <210> SEQ ID NO 342 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 342 tgtctagact ctgtggtatt gtgagga 27 <210> SEQ ID NO 343 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 343 tgtgtctgta aagactgaag ttgttgga 28 <210> SEQ ID NO 344 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 344 tgtgttttcc caagctcttc c 21 <210> SEQ ID NO 345 <211> LENGTH: 30 <212> TYPE: DNA

-6011	zinued	
<213> ORGANISM: Artificial Sequence <220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 345		
cgttccactt cccaagttta catcttcctc	30	
:210> SEQ ID NO 346		
<211> LENGTH: 25 <212> TYPE: DNA		
213> ORGANISM: Artificial Sequence		
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>		
:400> SEQUENCE: 346		
cgttccgcag actactatgg ctctc	25	
<210> SEQ ID NO 347 <211> LENGTH: 28		
<212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence <220> FEATURE:		
<223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 347		
tgttgaagtc aaaatacccc cagggggt	28	
<210> SEQ ID NO 348		
<211> LENGTH: 25 <212> TYPE: DNA		
<pre>&lt;212&gt; IIII. Dual </pre> <213> ORGANISM: Artificial Sequence		
<220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 348		
tgttgaggac ggctgctaat tgttg	25	
<210> SEQ ID NO 349		
<211> LENGTH: 31		
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence		
<220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 349		
~ cqttqtqata acaaaqtccc aatcatcqtt c	31	
<210> SEQ ID NO 350		
<211> LENGTH: 25 <212> TYPE: DNA		
<213> ORGANISM: Artificial Sequence		
<220> FEATURE: <223> OTHER INFORMATION: Primer		
<400> SEQUENCE: 350		
tgttgtgcac ttttggagaa tattt	25	
<210> SEQ ID NO 351		
<211> LENGTH: 25		
<212> TYPE: DNA <213> ORGANISM: Artificial Sequence		
<220> FEATURE:		
<223> OTHER INFORMATION: Primer		

	-continued
<pre>&lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (16)(16) &lt;223&gt; OTHER INFORMATION: I</pre>	
<400> SEQUENCE: 351	
tgttgttgat gagtetetge cagte	25
<210> SEQ ID NO 352 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 352	
tgttgttgat gagtctctgc cagtc	25
<210> SEQ ID NO 353 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 353	
tgtttatcca cggtcccact c	21
<pre>&lt;210&gt; SEQ ID NO 354 &lt;211&gt; LENGTH: 21 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;400&gt; SEQUENCE: 354</pre>	
tgtttatcca tggtcccact c	21
<pre>&lt;210&gt; SEQ ID NO 355 &lt;211&gt; LENGTH: 19 &lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer &lt;220&gt; FEATURE: &lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (11)(11) &lt;223&gt; OTHER INFORMATION: I</pre>	
<400> SEQUENCE: 355	
tgttttccca tgccctccc	19
<210> SEQ ID NO 356 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 356	
ttaggtcctg cacagccgca taatg	25
<210> SEQ ID NO 357 <211> LENGTH: 23 <212> TYPE: DNA	

	-continued
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 357	
tteteageta acaatttete age	23
<210> SEQ ID NO 358	
<211> LENGTH: 27 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 358	
ttctgcaatt acttgacacg acctcat	27
<210> SEQ ID NO 359	
<211> LENGTH: 25	
<pre>&lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence</pre>	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 359	
ttgcagagat ggagatcatt gtcca	25
<210> SEQ ID NO 360	
<211> LENGTH: 24 <212> TYPE: DNA	
<pre>&lt;213&gt; ORGANISM: Artificial Sequence &lt;220&gt; FEATURE:</pre>	
<2223> FLATORE: <223> OTHER INFORMATION: Primer	
<400> SEQUENCE: 360	
ttgcttggtt atcaccctga acca	24
<210> SEQ ID NO 361	
<211> LENGTH: 26	
<pre>&lt;212&gt; TYPE: DNA &lt;213&gt; ORGANISM: Artificial Sequence</pre>	
<pre>&lt;220&gt; FEATURE: &lt;223&gt; OTHER INFORMATION: Primer</pre>	
<400> SEQUENCE: 361	
ttgcttttta aaaatgcagc tgcatt	26
<210> SEQ ID NO 362	
<211> LENGTH: 26 <212> TYPE: DNA	
<213> ORGANISM: Artificial Sequence	
<220> FEATURE: <223> OTHER INFORMATION: Primer	
<220> FEATURE:	
<pre>&lt;221&gt; NAME/KEY: modified_base &lt;222&gt; LOCATION: (20)(21)</pre>	
<223> OTHER INFORMATION: I	
<400> SEQUENCE: 362	
ttgcttttta aaaatgcagc tgcatt	26
<210> SEQ ID NO 363	
<211> LENGTH: 26 <212> TYPE: DNA	
NATES TIPE: DNA	

-co			

<213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (20)...(21) <223> OTHER INFORMATION: I <400> SEQUENCE: 363 ttgcttttta aaaatgcagc tgcatt 26 <210> SEQ ID NO 364 <211> LENGTH: 26 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <220> FEATURE: <221> NAME/KEY: modified\_base <222> LOCATION: (18)...(21) <223> OTHER INFORMATION: I <400> SEQUENCE: 364 ttgcttttta aaaatgcagc tgcatt 26 <210> SEQ ID NO 365 <211> LENGTH: 24 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 365 ttggcggggt gaggaccttg aggg 24 <210> SEQ ID NO 366 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 366 ttgtcaatat catccagttg gc 22 <210> SEQ ID NO 367 <211> LENGTH: 21 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 367 ttgtgcacca tgcagttcat c 21 <210> SEQ ID NO 368 <211> LENGTH: 28 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 368 28 tttctgctcg tttataatgt ctacacat

-co	nti	inu	ιed

<210> SEO ID NO 369 <211> LENGTH: 23 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEOUENCE: 369 tttgagttgt gggcagtccc ttt 23 <210> SEQ ID NO 370 <211> LENGTH: 27 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 370 tttggtcctt gtcttatgtc cagaatg 27 <210> SEQ ID NO 371 <211> LENGTH: 22 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEQUENCE: 371 ttttggcgaa tattttgttt gg 22 <210> SEQ ID NO 372 <211> LENGTH: 25 <212> TYPE: DNA <213> ORGANISM: Artificial Sequence <220> FEATURE: <223> OTHER INFORMATION: Primer <400> SEOUENCE: 372 25 tgttccgcag accactatgg ctctc

What is claimed is:

1. An oligonucleotide primer 14 to 35 nucleobases in length comprising at least 70% sequence identity with SEQ ID NO: 47.

**2**. An oligonucleotide primer 14 to 35 nucleobases in length comprising at least 70% sequence identity with SEQ ID NO: 286.

3. A composition comprising the primer of claim 1.

**4**. The composition of claim 3 further comprising an oligonucleotide primer 14 to 35 nucleobases in length comprising at least 70% sequence identity with SEQ ID NO: 286.

**5**. The composition of claim 4 wherein either or both of said primers comprises at least one modified nucleobase.

**6**. The composition of claim 5 wherein said modified nucleobase is 5-propynyluracil or 5-propynylcytosine.

7. The composition of claim 4 wherein either or both of said primers comprises at least one universal nucleobase.

**8**. The composition of claim 7 wherein said universal nucleobase is inosine.

**9**. The composition of claim 4 wherein either or both of said primers further comprises a non-templated T residue on the 5'-end.

**10**. The composition of claim 4 wherein either or both of said primers comprises at least one non-template tag.

11. The composition of claim 4 wherein either or both of said primers comprises at least one molecular mass modi-fying tag.

12. A kit comprising the composition of claim 4.

**13**. The kit of claim 12 further comprising one or more primer pairs wherein each member of said one or more primer pairs is of a length of 14 to 35 nucleobases and has 70% to 100% sequence identity with the corresponding member from the group of primer pairs represented by SEQ ID NOs: 70:286, 165:286, 122:275, 100:336, and 61:324.

**14**. The kit of claim 12 further comprising at least one calibration polynucleotide.

**15**. The kit of claim 12 further comprising at least one anion exchange functional group linked to a magnetic bead.

**16**. A method for identification of a polyomavirus in a sample comprising:

- amplifying nucleic acid from said polyomavirus using the composition of claim 4 to obtain an amplification product;
- determining the molecular mass of said amplification product;
- optionally, determining the base composition of said amplification product from said molecular mass; and
- comparing said molecular mass or base composition with a plurality of molecular masses or base compositions of known polyomavirus identifying amplicons, wherein a match between said molecular mass or base composition and a member of said plurality of molecular masses or base compositions identifies said polyomavirus.

**17**. The method of claim 16 wherein said sample is a biological product.

**18**. A method of determining the presence or absence of a polyomavirus in a sample comprising:

- amplifying nucleic acid from said sample using the composition of claim 4 to obtain an amplification product;
- determining the molecular mass of said amplification product;
- optionally, determining the base composition of said amplification product from said molecular mass; and
- comparing said molecular mass or base composition of said amplification product with the known molecular masses or base compositions of one or more known

polyomavirus identifying amplicons, wherein a match between said molecular mass or base composition of said amplification product and the molecular mass or base composition of one or more known polyomavirus identifying amplicons indicates the presence of said polyomavirus in said sample.

**19**. The method of claim 18 wherein said sample comprises a biological product.

**20**. A method for determination of the quantity of an unknown polyomavirus in a sample comprising:

- contacting said sample with the composition of claim 4 and a known quantity of a calibration polynucleotide comprising a calibration sequence;
- concurrently amplifying nucleic acid from said unknown polyomavirus and nucleic acid from said calibration polynucleotide in said sample with the composition of claim 4 to obtain a first amplification product comprising a polyomavirus identifying amplicon and a second amplification product comprising a calibration amplicon;
- determining the molecular mass and abundance for said polyomavirus identifying amplicon and said calibration amplicon; and
- distinguishing said polyomavirus identifying amplicon from said calibration amplicon based on molecular mass, wherein comparison of polyomavirus identifying amplicon abundance and calibration amplicon abundance indicates the quantity of polyomavirus in said sample.

\* \* \* \* \*