

(19) **DANMARK**

(10) **DK/EP 3250592 T3**



(12) **Oversættelse af  
europæisk patentskrift**

Patent- og  
Varemærkestyrelsen

- 
- (51) Int.Cl.: **C 07 K 16/18 (2006.01)**
- (45) Oversættelsen bekendtgjort den: **2023-01-09**
- (80) Dato for Den Europæiske Patentmyndigheds bekendtgørelse om meddelelse af patentet: **2022-10-05**
- (86) Europæisk ansøgning nr.: **16702812.5**
- (86) Europæisk indleveringsdag: **2016-01-28**
- (87) Den europæiske ansøgnings publiceringsdag: **2017-12-06**
- (86) International ansøgning nr.: **IB2016050415**
- (87) Internationalt publikationsnr.: **WO2016120810**
- (30) Prioritet: **2015-01-28 US 201562109002 P** **2015-12-11 US 201562266556 P**
- (84) Designerede stater: **AL AT BE BG CH CY CZ DE DK EE ES FI FR GB GR HR HU IE IS IT LI LT LU LV MC MK MT NL NO PL PT RO RS SE SI SK SM TR**
- (73) Patenthaver: **Prothema Biosciences Limited, 77 Sir John Rogerson's Quay, Block C , Grand Canal Docklands, Dublin 2, D02 VK60, Irland**  
**University Health Network, 190 Elizabeth Street , R. Fraser Elliott Building - Room 1S-417, Toronto, Ontario M5G 2C4, Canada**
- (72) Opfinder: **LIU, Yue, 815 Peary Lane, Foster City, California 94404, USA**  
**NIJJAR, Tarlochan S., 712 Ironbark Court, Orinda, California 94563, USA**  
**CHAKRABARTTY, Avijit, 77 Points O'Woods Drive, Vaughn, Ontario L4K 2E1, Canada**  
**HIGAKI, Jeffrey, N., 2071 Queens Lane, San Mateo, CA 94402, USA**
- (74) Fuldmægtig i Danmark: **Budde Schou A/S, Dronningens Tværgade 30, 1302 København K, Danmark**
- (54) Benævnelse: **ANTI-TRANSTHYRETIN-ANTISTOFFER**
- (56) Fremdragne publikationer:  
**WO-A2-2014/124334**  
**WO-A2-2015/010118**  
**US-A1- 2014 056 904**  
**MONICHAN PHAY ET AL: "Transthyretin Aggregate-Specific Antibodies Recognize Cryptic Epitopes on Patient-Derived Amyloid Fibrils", REJUVENATION RESEARCH, vol. 17, no. 2, 1 April 2014 (2014-04-01), pages 97-104, XP055256411, US ISSN: 1549-1684, DOI: 10.1089/rej.2013.1524**  
**GOLDSTEINS G ET AL: "Exposure of cryptic epitopes on transthyretin only in amyloid and in amyloidogenic mutants", PROCEEDINGS OF THE NATIONAL ACADEMY OF SCIENCES, NATIONAL ACADEMY OF SCIENCES, US, vol. 96, no. 6, 16 March 1999 (1999-03-16) , pages 3108-3113, XP002987251, ISSN: 0027-8424, DOI: 10.1073/PNAS.96.6.3108**  
**Olivier Léger ET AL: "Antibody Drug Discovery Chapter 1: "Humanization of Antibodies"" In: "Molecular Medicine and Medicinal Chemistry", 1 January 2011 (2011-01-01), XP055119233, pages 1-23, the whole document**

Fortsættes ...

**ALMAGRO JUAN C ET AL: "Humanization of antibodies", FRONTIERS IN BIOSCIENCE, FRONTIERS IN BIOSCIENCE, ALBERTSON, NY, US, vol. 13, 1 January 2008 (2008-01-01), pages 1619-1633, XP009126790, ISSN: 1093-9946**

**Description****REFERENCE TO A SEQUENCE LISTING**

5 **[0001]** This application includes an electronic sequence listing in a file named 473373\_SEQLIST.TXT, created on January 28, 2016, and containing 135,083 bytes, which is hereby incorporated by reference in its entirety for all purposes.

**BACKGROUND**

10 **[0002]** Several diseases are thought to be caused by the abnormal folding and aggregation of disease-specific proteins. These proteins can accumulate into pathologically diagnostic accumulations, known as amyloids, which are visualized by certain histologic stains. Amyloids are thought to elicit inflammatory responses and have multiple negative consequences for the involved tissues. In addition, smaller aggregates of abnormally folded protein may exist and exert cytotoxic effects.

15 **[0003]** Transthyretin (TTR) is one of the many proteins that are known to misfold and aggregate (e.g., undergo amyloidogenesis). Transthyretin-related amyloidosis encompasses two forms of disease: familial disease arising from misfolding of a mutated or variant TTR, and a sporadic, non-genetic disease caused by misaggregation of wild-type TTR. The process of TTR amyloidogenesis can cause pathology in the nervous system and/or heart, as well as in other tissues.

20 **[0004]** US2014056904 discusses antibodies binding to misfolded TTR, including a polyclonal antibody raised against residues 89-97 of TTR, and diagnostic and therapeutic applications.

**[0005]** WO2014124334 discusses antibodies selectively binding to non-native TTR and diagnostic uses thereof.

**[0006]** WO2015010118 discusses isolation of human antibodies with catalytic activity in dissolving misfolded TTR.

**[0007]** Phay et al. (Rejuvenation Research, 2014, 17, 97-104) discusses antibodies binding to transthyretin amyloid deposits in cardiac tissue sections from two familial amyloidotic polyneuropathy patients.

25 **[0008]** Goldsteins et al. (PNAS, 1999, 96, 3108-3113) discusses formation of several monoclonal antibodies against a TTR mutant with amyloid-forming potential. The antibodies were reported to be unable to inhibit amyloid formation.

**SUMMARY**

30 **[0009]** In one aspect the present invention provides an antibody that binds to transthyretin and comprises a mature heavy chain variable region of SEQ ID NO: 11 and a mature light chain variable region of SEQ ID NO: 19, or a mature heavy chain variable region of SEQ ID NO: 65 and a mature light chain variable region of SEQ ID NO:76, as set out in the claims.

35 **[0010]** The present invention also provides for a pharmaceutical composition comprising such an antibody, and a nucleic acid or nucleic acids encoding the heavy chain and light chain of such an antibody, as set out in the claims.

**[0011]** A further aspect of the present invention is a recombinant expression vector comprising such a nucleic acid as set out in the claims, and a host cell transformed with said recombinant expression vector as set out in the claims.

40 **[0012]** In another aspect, the present invention provides a method of producing an antibody as set out in the claims the method comprising: culturing cells transformed with nucleic acids encoding the heavy and light chains of the antibody, so that the cells secrete the antibody; and purifying the antibody from cell culture media.

**[0013]** In yet a further aspect the present invention provides an antibody as set out in the claims, for use in treating or effecting prophylaxis of a transthyretin-mediated amyloidosis, or delaying the onset of a transthyretin-mediated amyloidosis in the subject.

45 **[0014]** In a still further aspect, the present invention provides an in vitro method of diagnosing a transthyretin-mediated amyloidosis in a subject, comprising contacting a biological sample from the subject with an effective amount of an antibody as set out in the claims.

50 **[0015]** Disclosed herein are antibodies that bind to transthyretin and comprise three heavy chain CDRs and three light chain CDRs substantially from antibody 14G8. Some antibodies comprise three Kabat heavy chain CDRs and three Kabat light chain CDRs of the antibody 14G8, except that positions H52 and L26 can each be N or S. Some antibodies comprise three Kabat heavy chain CDRs and three Kabat light chain CDRs of the antibody 14G8. Optionally, CDR-H1 is a composite Kabat-Chothia CDR of the antibody 14G8. Some antibodies comprise three Kabat heavy chain CDRs of the antibody 9D5. Optionally, CDR-H1 is a composite Kabat-Chothia CDR of the antibody 9D5.

55 **[0016]** Some antibodies bind to the same epitope on transthyretin as 9D5 or 14G8 and comprises three light chain CDRs and three heavy chain CDRs, wherein: (a) each CDR has at least 90% sequence identity to a corresponding CDR from the heavy and light chain variable regions of 9D5; or (b) each CDR has at least 90% sequence identity to a corresponding CDR from the heavy and light chain variable regions of 14G8, except that position L26 can be N or S. Some antibodies comprise: (a) the three heavy chain CDRs and the three light chain CDRs of 9D5; or (b) the three heavy chain CDRs and the three light chain CDRs of 14G8, except that position L26 can be N or S. In some antibodies,

the three heavy chain CDRs and the three light chain CDRs of 9D5 are SEQ ID NOS: 13-15 and 24-26, respectively. In some antibodies, the three heavy chain CDRs and the three light chain CDRs of 14G8 are SEQ ID NOS: 67-69 and 77-79, respectively, except that position L26 can be N or S.

**[0017]** Any of the above antibodies can be a monoclonal antibody. Any of the above antibodies can be a chimeric, humanized, veneered, or human antibody. Any of the above antibodies can have human IgG1 isotype, human IgG2 isotype, or human IgG4 isotype.

**[0018]** Another aspect of the disclosure provides humanized or chimeric antibodies of mouse antibody that specifically binds to transthyretin, wherein the mouse antibody is characterized by a mature heavy chain variable region of SEQ ID NO: 61 and a mature light chain variable region of SEQ ID NO: 70, except that position H52 can be S or N, position H69 can be F or I, position L26 can be S or N, and the R at position L107 is optional. Optionally, the antibodies are a humanized or chimeric 9D5 or 14G8 antibody that specifically binds to transthyretin, wherein 9D5 is a mouse antibody characterized by a mature heavy chain variable region of SEQ ID NO: 1 and a mature light chain variable region of SEQ ID NO: 16, and 14G8 is a mouse antibody characterized by a mature heavy chain variable region of SEQ ID NO: 61 and a mature light chain variable region of SEQ ID NO: 70. Optionally, the humanized antibodies comprise: (a) a humanized mature heavy chain variable region comprising the three heavy chain CDRs of 9D5 and a humanized mature light chain variable region comprising the three light chain CDRs of 9D5; or (b) a humanized mature heavy chain variable region comprising the three heavy chain CDRs of 14G8 and a humanized mature light chain variable region comprising the three light chain CDRs of 14G8, except that position L26 can be N or S. Optionally, the three heavy chain CDRs and the three light chain CDRs of 9D5 are SEQ ID NOS: 13-15 and 24-26, respectively. Optionally, the three heavy chain CDRs and the three light chain CDRs of 14G8 are SEQ ID NOS: 67-69 and 77-79, respectively, except that position L26 can be N or S. Optionally, any differences in CDRs of the mature heavy chain variable region and mature light chain variable region from SEQ ID NOS: 1 and 16, respectively, reside in positions H60-H65. Optionally, any differences in CDRs of the mature heavy chain variable region and mature light chain variable region from SEQ ID NOS: 61 and 70, respectively, reside in positions H60-H65, except that position L26 can be N or S.

**[0019]** Optionally, the humanized antibody comprises a humanized mature heavy chain variable region having an amino acid sequence at least 90% identical to any one of SEQ ID NOS: 5-12 and a humanized mature light chain variable region having an amino acid sequence at least 90% identical to any one of SEQ ID NOS: 19-23, except that position H19 can be R or K, position H40 can be A or T, position H44 can be G or R, position H49 can be S or A, position H77 can be S or T, position H82a can be N or S, position H83 can be R or K, position H84 can be A or S, and position H89 can be V or M. Optionally, at least one of the following positions is occupied by the amino acid as specified: position H42 is occupied by E, position H47 is occupied by L, position H69 is occupied by F, position H82 is occupied by S, position H82b is occupied by L, position H108 is occupied by L, position L8 is occupied by A, position L9 is occupied by P, position L18 is occupied by S, position L19 is occupied by V, position L36 is occupied by F, position L39 is occupied by R, position L60 is occupied by S, position L70 is occupied by A, and position L74 is occupied by R. Optionally, positions H47, H69, and H82 are occupied by L, F, and S, respectively. Optionally, positions H47, H69, H82, and H82b are occupied by L, F, S, and L, respectively. Optionally, positions H42, H47, and H108 are occupied by E, L, and L, respectively. Optionally, positions H69, H82, and H82b are occupied by F, S, and L, respectively. Optionally, positions H47 and H108 are each occupied by L. Optionally, positions H82 and H82b are occupied by S and L, respectively. Optionally, positions H42, H47, and H82b are occupied by E, L, and L, respectively. Optionally, position L36 is occupied by F. Optionally, position L60 is occupied by S. Optionally, positions L8, L9, L19, L36, L39, L60, L70, and L74 are occupied by A, P, V, F, R, S, A, and R, respectively. Optionally, positions L8, L9, L18, L19, L36, L39, L60, L70, and L74 are occupied by A, P, S, V, F, R, S, A, and R, respectively.

**[0020]** Optionally, the humanized antibody comprises a mature heavy chain variable region having an amino acid sequence at least 95% identical to any one of SEQ ID NOS: 5-12 and a mature light chain variable region having an amino acid sequence at least 95% identical to any one of SEQ ID NOS: 19-23, except that position H19 can be R or K, position H40 can be A or T, position H44 can be G or R, position H49 can be S or A, position H77 can be S or T, position H82a can be N or S, position H83 can be R or K, position H84 can be A or S, and position H89 can be V or M. Optionally, the humanized antibody comprises a mature heavy chain variable region having an amino acid sequence at least 98% identical to any one of SEQ ID NOS: 5-12 and a mature light chain variable region having an amino acid sequence at least 98% identical to any one of SEQ ID NOS: 19-23, except that position H19 can be R or K, position H40 can be A or T, position H44 can be G or R, position H49 can be S or A, position H77 can be S or T, position H82a can be N or S, position H83 can be R or K, position H84 can be A or S, and position H89 can be V or M. Optionally, the mature heavy chain variable region has an amino acid sequence of any one of SEQ ID NO: 5-12 and the mature light chain variable region has an amino acid sequence of any one of SEQ ID NO: 19-23. In some antibodies, the mature heavy chain variable region has an amino acid sequence of SEQ ID NO: 11 and the mature light chain variable region has an amino acid sequence of SEQ ID NO: 19.

**[0021]** Optionally, the humanized antibody comprises a humanized mature heavy chain variable region having an amino acid sequence at least 90% identical to any one of SEQ ID NOS: 64-66 and a humanized mature light chain

variable region having an amino acid sequence at least 90% identical to any one of SEQ ID NOS: 74-76, except that position H82a can be N or S, position H83 can be R or K, position H84 can be A or S, position H89 can be V or M, and position L18 can be S or P. Optionally, at least one of the following positions is occupied by the amino acid as specified: position H1 is occupied by E, position H47 is occupied by L, and position L36 is occupied by F. Optionally, positions H1 and H47 are occupied by E and L, respectively. Optionally, position L36 is occupied by F. Optionally, at least one of the following positions is occupied by the amino acid as specified: position H3 is occupied by K, position H105 is occupied by T, position L8 is occupied by A, position L9 is occupied by P, position L19 is occupied by V, position L26 is occupied by S, position L60 is occupied by S, and position L70 is occupied by A. Optionally, positions H3 and H105 are occupied by K and T, respectively. Optionally, positions L8, L9, L19, and L70 are occupied by A, P, V, and A, respectively. Optionally, positions L26 and L60 are each occupied by S.

**[0022]** Optionally, the humanized antibody comprises a mature heavy chain variable region having an amino acid sequence at least 95% identical to any one of SEQ ID NOS: 64-66 and a mature light chain variable region having an amino acid sequence at least 95% identical to any one of SEQ ID NOS: 74-76, except that position H82a can be N or S, position H83 can be R or K, position H84 can be A or S, position H89 can be V or M, and position L18 can be S or P. Optionally, the humanized antibody comprises a mature heavy chain variable region having an amino acid sequence at least 98% identical to any one of SEQ ID NOS: 64-66 and a mature light chain variable region having an amino acid sequence at least 98% identical to any one of SEQ ID NOS: 74-76, except that position H82a can be N or S, position H83 can be R or K, position H84 can be A or S, position H89 can be V or M, and position L18 can be S or P. Optionally, the mature heavy chain variable region has an amino acid sequence of any one of SEQ ID NOS: 64-66, and the mature light chain variable region has an amino acid sequence of any one of SEQ ID NOS: 74-76. In some antibodies, the mature heavy chain variable region has an amino acid sequence of SEQ ID NO: 65 and the mature light chain variable region has an amino acid sequence of SEQ ID NO: 76.

**[0023]** Any of the above antibodies can be an intact antibody, a binding fragment, a single-chain antibody, a Fab, or a Fab'2 fragment. In any of the above antibodies, the mature light chain variable region can be fused to a light chain constant region and the mature heavy chain variable region can be fused to a heavy chain constant region. Optionally, the heavy chain constant region is a mutant form of a natural human heavy chain constant region which has reduced binding to a Fc $\gamma$  receptor relative to the natural human heavy chain constant region. Optionally, the heavy chain constant region is of IgG1 isotype. Optionally, the mature heavy chain variable region is fused to a heavy chain constant region having the sequence of SEQ ID NO: 103 and/or the mature light chain variable region is fused to a light chain constant region having the sequence of SEQ ID NO: 104 or 105.

**[0024]** Also disclosed herein are pharmaceutical compositions comprising any of the antibodies and a pharmaceutically acceptable carrier.

**[0025]** In a further aspect of the disclosure, provided herein are nucleic acids encoding the heavy chain and/or light chain of any of the above antibodies, such as any one of SEQ ID NOS: 40, 42, 44-56, 87, 89, 91-96, and 106-108.

**[0026]** Also disclosed is a recombinant expression vector comprising a nucleic acid as described above, and a host cell transformed with the recombinant expression vector.

**[0027]** The disclosure further provides a method of humanizing an antibody, the method comprising: (a) selecting one or more acceptor antibodies; (b) identifying the amino acid residues of the mouse antibody to be retained; (c) synthesizing a nucleic acid encoding a humanized heavy chain comprising CDRs of the mouse antibody heavy chain and a nucleic acid encoding a humanized light chain comprising CDRs of the mouse antibody light chain; and (d) expressing the nucleic acids in a host cell to produce a humanized antibody; wherein the mouse antibody is 9D5 or 14G8, wherein 9D5 is characterized by a mature heavy chain variable region of SEQ ID NO: 1 and a mature light chain variable region of SEQ ID NO: 16, and 14G8 is characterized by a mature heavy chain variable region of SEQ ID NO: 61 and a mature light chain variable region of SEQ ID NO: 70.

**[0028]** The disclosure further provides a method of producing a humanized, chimeric, or veneered antibody, the method comprising: (a) culturing cells transformed with nucleic acids encoding the heavy and light chains of the antibody, so that the cells secrete the antibody; and (b) purifying the antibody from cell culture media; wherein the antibody is a humanized, chimeric, or veneered form of 9D5 or 14G8.

**[0029]** The disclosure further provides a method of producing a cell line producing a humanized, chimeric, or veneered antibody, the method comprising: (a) introducing a vector encoding heavy and light chains of an antibody and a selectable marker into cells; (b) propagating the cells under conditions to select for cells having increased copy number of the vector; (c) isolating single cells from the selected cells; and (d) banking cells cloned from a single cell selected based on yield of antibody; wherein the antibody is a humanized, chimeric, or veneered form of 9D5 or 14G8. Optionally, the method further comprises propagating the cells under selective conditions and screening for cell lines naturally expressing and secreting at least 100 mg/L/10<sup>6</sup> cells/24h.

**[0030]** The disclosure further provides a method of inhibiting or reducing aggregation of transthyretin in a subject having or at risk of developing a transthyretin-mediated amyloidosis, comprising administering to the subject an effective regime of any of the above antibodies, thereby inhibiting or reducing aggregation of transthyretin in the subject.

**[0031]** The disclosure further provides a method of inhibiting or reducing transthyretin fibril formation in a subject having or at risk of developing a transthyretin-mediated amyloidosis, comprising administering to the subject an effective regime of any of the above antibodies, thereby inhibiting or reducing transthyretin accumulation in the subject.

**[0032]** The disclosure further provides a method of reducing transthyretin deposits in a subject having or at risk of developing a transthyretin-mediated amyloidosis, comprising administering to the subject an effective regime of any of the above antibodies, thereby reducing transthyretin deposits in the subject.

**[0033]** The disclosure further provides a method of clearing aggregated transthyretin in a subject having or at risk of developing a transthyretin-mediated amyloidosis, comprising administering to the subject an effective regime of any of the above antibodies, thereby clearing aggregated transthyretin from the subject relative to a subject having or at risk of developing a transthyretin-mediated amyloidosis who has not received the antibody.

**[0034]** The disclosure further provides a method of stabilizing a non-toxic conformation of transthyretin in a subject having or at risk of developing a transthyretin-mediated amyloidosis, comprising administering to the subject an effective regime of any of the above antibodies, thereby stabilizing a non-toxic conformation of transthyretin in the subject.

**[0035]** The disclosure further provides a method of treating or effecting prophylaxis of a transthyretin-mediated amyloidosis in a subject, comprising administering to the subject an effective regime of any of the above antibodies.

**[0036]** The disclosure further provides a method of delaying the onset of a transthyretin-mediated amyloidosis in a subject, comprising administering to the subject an effective regime of any of the above antibodies.

**[0037]** The disclosure further provides a method of diagnosing a transthyretin-mediated amyloidosis in a subject, comprising contacting a biological sample from the subject with an effective amount of any of the above antibodies.

Optionally, the method further comprises detecting the binding of antibody to transthyretin, wherein the presence of bound antibody indicates the subject has a transthyretin-mediated amyloidosis. Optionally, the method further comprises comparing binding of the antibody to the biological sample with binding of the antibody to a control sample, whereby increased binding of the antibody to the biological sample relative to the control sample indicates the subject has a transthyretin-mediated amyloidosis. Optionally, the biological sample and the control sample comprise cells of the same tissue origin. Optionally, the biological sample and/or the control sample is blood, serum, plasma, or solid tissue. Optionally, the solid tissue is from the heart, peripheral nervous system, autonomic nervous system, kidneys, eyes, or gastrointestinal tract.

**[0038]** In any of the above methods, the transthyretin-mediated amyloidosis can optionally be a familial transthyretin amyloidosis or a sporadic transthyretin amyloidosis. Optionally, the familial transthyretin amyloidosis is familial amyloid cardiomyopathy (FAC), familial amyloid polyneuropathy (FAP), or central nervous system selective amyloidosis (CNSA). Optionally, the sporadic transthyretin amyloidosis is senile systemic amyloidosis (SSA) or senile cardiac amyloidosis (SCA). In any of the above methods, the transthyretin-mediated amyloidosis can optionally be associated with amyloid accumulation in the heart, peripheral nervous system, autonomic nervous system, kidneys, eyes, or gastrointestinal tract of the subject.

**[0039]** The disclosure further provides a method of detecting the presence or absence of transthyretin deposits in a subject, comprising contacting a biological sample from the subject suspected of comprising the amyloid accumulation with an effective amount of any of the above antibodies. Optionally, the method further comprises detecting the binding of antibody to transthyretin, wherein detection of bound antibody indicates the presence of transthyretin deposits. Optionally, the method further comprises comparing binding of the antibody to the biological sample with binding of the antibody to a control sample, whereby increased binding of the antibody to the biological sample relative to the control sample indicates the subject has a transthyretin-mediated amyloidosis. Optionally, the biological sample and the control sample comprise cells of the same tissue origin. Optionally, the biological sample and/or the control sample is blood, serum, plasma, or solid tissue. Optionally, the solid tissue is from the heart, peripheral nervous system, autonomic nervous system, kidneys, eyes, or gastrointestinal tract.

**[0040]** The disclosure further provides a method of determining a level of transthyretin deposits in a subject, comprising administering any of the above antibodies and detecting the presence of bound antibody in the subject. Optionally, the presence of bound antibody is determined by positron emission tomography (PET).

## BRIEF DESCRIPTION OF THE DRAWINGS

### [0041]

FIG. 1A depicts an alignment of heavy chain variable regions of the mouse 9D5 antibody, mouse model antibodies, human acceptor antibodies, and humanized versions of the 9D5 antibody. The CDRs as defined by Kabat are enclosed in boxes, except that the first enclosed box is a composite of the Chothia CDR-H1 and the Kabat CDR-H1, with the Kabat CDR-H1 underlined and bolded.

FIG. 1B depicts an alignment of light chain variable regions of the mouse 9D5 antibody, a mouse model antibody,

a human acceptor antibody, and humanized versions of the 9D5 antibody. The CDRs as defined by Kabat are enclosed in boxes.

5 FIG. 2A depicts an alignment of heavy chain variable regions of the mouse 14G8 antibody, a mouse model antibody, human acceptor antibodies, and humanized versions of the 14G8 antibody. The CDRs as defined by Kabat are enclosed in boxes, except that the first enclosed box is a composite of the Chothia CDR-H1 and the Kabat CDR-H1, with the Kabat CDR-H1 underlined and bolded.

10 FIG. 2B depicts an alignment of light chain variable regions of the mouse 14G8 antibody, a mouse model antibody, human acceptor antibodies, and humanized versions of the 14G8 antibody. The CDRs as defined by Kabat are enclosed in boxes.

15 FIG. 3A & 3B: FIG. 3A depicts binding curves of murine 5A1, 6C1, 9D5, and 14G8 antibodies to pH4-treated TTR. FIG. 3B depicts binding curves of murine 5A1, 6C1, 9D5, and 14G8 antibodies to native TTR.

FIG. 4A, 4B & 4C: FIG. 4A depicts the inhibition of TTR-Y78F fiber formation by mis-TTR antibodies. FIG. 4B depicts the inhibition of TTR-V122I fiber formation by 14G8. FIG. 4C depicts the inhibition of TTR-V122I fiber formation by a control antibody.

20 FIG. 5A & 5B: FIG. 5A depicts a densitometry analysis of a Western blot analysis of plasma samples from patients confirmed for V30M ATTR (Sample #21, #22, #23, #24, #25, #27) and samples from normal subjects (Sample #11, #12, #15, #18, #19, #20) using the 9D5 mis-TTR antibody. FIG. 5B depicts a densitometry analysis of a Western blot analysis of the same samples using the 5A1 mis-TTR antibody.

25 FIG. 6 depicts a MesoScale Discovery (MSD) plate assay of plasma samples from patients confirmed for V30M ATTR (#21, #22, #23, #24, #25, #27) and samples from normal subjects (Sample #11, #12, #15, #18, #19, #20) using the 6C1 antibody.

30 FIGS. 7A & 7B: FIG. 7A depicts the effect of antibody 14G8 on the uptake of F87M/L110M TTR by THP-1 cells. FIG. 7B depicts the effect of each of the mis-TTR antibodies on the uptake of V30M TTR by THP-1 cells.

## BRIEF DESCRIPTION OF THE SEQUENCES

### [0042]

35 SEQ ID NO: 1 sets forth the amino acid sequence of the heavy chain variable region of the mouse 9D5 antibody.

SEQ ID NO: 2 sets forth the amino acid sequence of the mouse heavy chain variable region structure template 1SEQ H.

40 SEQ ID NO: 3 sets forth the amino acid sequence of the heavy chain variable acceptor ACC# BAC02114.

SEQ ID NO: 4 sets forth the amino acid sequence of the heavy chain variable acceptor ACC#AAX82494.1.

45 SEQ ID NO: 5 sets forth the amino acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 1 (Hu9D5VHv1).

50 SEQ ID NO: 6 sets forth the amino acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 2 (Hu9D5VHv2).

SEQ ID NO: 7 sets forth the amino acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 2b (Hu9D5VHv2b).

55 SEQ ID NO: 8 sets forth the amino acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 3 (Hu9D5VHv3).

SEQ ID NO: 9 sets forth the amino acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 3b (Hu9D5VHv3b).

SEQ ID NO: 10 sets forth the amino acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 4 (Hu9D5VHv4).

5 SEQ ID NO: 11 sets forth the amino acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 4b (Hu9D5VHv4b).

SEQ ID NO: 12 sets forth the amino acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 5 (Hu9D5VHv5).

10 SEQ ID NO: 13 sets forth the amino acid sequence of Kabat CDR-H1 of the mouse 9D5 antibody.

SEQ ID NO: 14 sets forth the amino acid sequence of Kabat CDR-H2 of the mouse 9D5 antibody.

15 SEQ ID NO: 15 sets forth the amino acid sequence of Kabat CDR-H3 of the mouse 9D5 antibody.

SEQ ID NO: 16 sets forth the amino acid sequence of the light chain variable region of the mouse 9D5 antibody.

20 SEQ ID NO: 17 sets forth the amino acid sequence of the mouse light chain variable region structure template 1MJU\_L.

SEQ ID NO: 18 sets forth the amino acid sequence of the light chain variable acceptor ACC# ABC66952.

25 SEQ ID NO: 19 sets forth the amino acid sequence of the light chain variable region of the humanized 9D5 antibody version 1 (Hu9D5VLv1).

SEQ ID NO: 20 sets forth the amino acid sequence of the light chain variable region of the humanized 9D5 antibody version 2 (Hu9D5VLv2).

30 SEQ ID NO: 21 sets forth the amino acid sequence of the light chain variable region of the humanized 9D5 antibody version 3 (Hu9D5VLv3).

SEQ ID NO: 22 sets forth the amino acid sequence of the light chain variable region of the humanized 9D5 antibody version 4 (Hu9D5VLv4).

35 SEQ ID NO: 23 sets forth the amino acid sequence of the light chain variable region of the humanized 9D5 antibody version 5 (Hu9D5VLv5).

SEQ ID NO: 24 sets forth the amino acid sequence of Kabat CDR-L1 of the mouse 9D5 antibody.

40 SEQ ID NO: 25 sets forth the amino acid sequence of Kabat CDR-L2 of the mouse 9D5 antibody.

SEQ ID NO: 26 sets forth the amino acid sequence of Kabat CDR-L3 of the mouse 9D5 antibody.

45 SEQ ID NO: 27 sets forth the amino acid sequence of humanized 9D5 heavy chain version 1.

SEQ ID NO: 28 sets forth the amino acid sequence of humanized 9D5 heavy chain version 2.

SEQ ID NO: 29 sets forth the amino acid sequence of humanized 9D5 heavy chain version 2b.

50 SEQ ID NO: 30 sets forth the amino acid sequence of humanized 9D5 heavy chain version 3.

SEQ ID NO: 31 sets forth the amino acid sequence of humanized 9D5 heavy chain version 3b.

55 SEQ ID NO: 32 sets forth the amino acid sequence of humanized 9D5 heavy chain version 4.

SEQ ID NO: 33 sets forth the amino acid sequence of humanized 9D5 heavy chain version 4b.

SEQ ID NO: 34 sets forth the amino acid sequence humanized 9D5 heavy chain version 5.

- SEQ ID NO: 35 sets forth the amino acid sequence of humanized 9D5 light chain version 1.
- SEQ ID NO: 36 sets forth the amino acid sequence of humanized 9D5 light chain version 2.
- 5 SEQ ID NO: 37 sets forth the amino acid sequence of humanized 9D5 light chain version 3.
- SEQ ID NO: 38 sets forth the amino acid sequence of humanized 9D5 light chain version 4.
- 10 SEQ ID NO: 39 sets forth the amino acid sequence of humanized 9D5 light chain version 5.
- SEQ ID NO: 40 sets forth the nucleic acid sequence of the heavy chain variable region of the mouse 9D5 antibody with signal peptide.
- 15 SEQ ID NO: 41 sets forth the amino acid sequence of the heavy chain variable region of the mouse 9D5 antibody with signal peptide.
- SEQ ID NO: 42 sets forth the nucleic acid sequence of the light chain variable region of the mouse 9D5 antibody with signal peptide.
- 20 SEQ ID NO: 43 sets forth the amino acid sequence of the light chain variable region of the mouse 9D5 antibody with signal peptide.
- SEQ ID NO: 44 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 1 (Hu9D5VHv1).
- 25 SEQ ID NO: 45 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 2 (Hu9D5VHv2).
- SEQ ID NO: 46 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 2b (Hu9D5VHv2b).
- 30 SEQ ID NO: 47 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 3 (Hu9D5VHv3).
- SEQ ID NO: 48 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 3b (Hu9D5VHv3b).
- 35 SEQ ID NO: 49 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 4 (Hu9D5VHv4).
- 40 SEQ ID NO: 50 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 4b (Hu9D5VHv4b).
- SEQ ID NO: 51 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 9D5 antibody version 5 (Hu9D5VHv5).
- 45 SEQ ID NO: 52 sets forth the nucleic acid sequence of the light chain variable region of the humanized 9D5 antibody version 1 (Hu9D5VLv1).
- SEQ ID NO: 53 sets forth the nucleic acid sequence of the light chain variable region of the humanized 9D5 antibody version 2 (Hu9D5VLv2).
- 50 SEQ ID NO: 54 sets forth the nucleic acid sequence of the light chain variable region of the humanized 9D5 antibody version 3 (Hu9D5VLv3).
- 55 SEQ ID NO: 55 sets forth the nucleic acid sequence of the light chain variable region of the humanized 9D5 antibody version 4 (Hu9D5VLv4).

SEQ ID NO: 56 sets forth the nucleic acid sequence of the light chain variable region of the humanized 9D5 antibody version 5 (Hu9D5VLv5).

5 SEQ ID NO: 57 sets forth the amino acid sequence of the mouse 9D5 heavy chain variable region signal peptide.

SEQ ID NO: 58 sets forth the nucleic acid sequence of the mouse 9D5 heavy chain variable region signal peptide.

SEQ ID NO: 59 sets forth the amino acid sequence of the mouse 9D5 light chain variable region signal peptide.

10 SEQ ID NO: 60 sets forth the nucleic acid sequence of the mouse 9D5 light chain variable region signal peptide.

SEQ ID NO: 61 sets forth the amino acid sequence of the heavy chain variable region of the mouse 14G8 antibody.

15 SEQ ID NO: 62 sets forth the amino acid sequence of the mouse heavy chain variable region structure template 1MQK\_H.

SEQ ID NO: 63 sets forth the amino acid sequence of the heavy chain variable acceptor ACC# AAD30410.1.

20 SEQ ID NO: 64 sets forth the amino acid sequence of the heavy chain variable region of the humanized 14G8 antibody version 1 (Hu14G8VHv1).

SEQ ID NO: 65 sets forth the amino acid sequence of the heavy chain variable region of the humanized 14G8 antibody version 2 (Hu14G8VHv2).

25 SEQ ID NO: 66 sets forth the amino acid sequence of the heavy chain variable region of the humanized 14G8 antibody version 3 (Hu14G8VHv3).

SEQ ID NO: 67 sets forth the amino acid sequence of Kabat CDR-H1 of the mouse 14G8 antibody.

30 SEQ ID NO: 68 sets forth the amino acid sequence of Kabat CDR-H2 of the mouse 14G8 antibody.

SEQ ID NO: 69 sets forth the amino acid sequence of Kabat CDR-H3 of the mouse 14G8 antibody.

35 SEQ ID NO: 70 sets forth the amino acid sequence of the light chain variable region of the mouse 14G8 antibody.

SEQ ID NO: 71 sets forth the amino acid sequence of the mouse light chain variable region structure template 1MJU\_L.

40 SEQ ID NO: 72 sets forth the amino acid sequence of the light chain variable acceptor ACC# ABA71374.1.

SEQ ID NO: 73 sets forth the amino acid sequence of the light chain variable acceptor ACC# ABC66952.1.

45 SEQ ID NO: 74 sets forth the amino acid sequence of the light chain variable region of the humanized 14G8 antibody version 1 (Hu14G8VLv1).

SEQ ID NO: 75 sets forth the amino acid sequence of the light chain variable region of the humanized 14G8 antibody version 2 (Hu14G8VLv2).

50 SEQ ID NO: 76 sets forth the amino acid sequence of the light chain variable region of the humanized 14G8 antibody version 3 (Hu14G8VLv3).

SEQ ID NO: 77 sets forth the amino acid sequence of Kabat CDR-L1 of the mouse 14G8 antibody.

55 SEQ ID NO: 78 sets forth the amino acid sequence of Kabat CDR-L2 of the mouse 14G8 antibody.

SEQ ID NO: 79 sets forth the amino acid sequence of Kabat CDR-L3 of the mouse 14G8 antibody.

SEQ ID NO: 80 sets forth the amino acid sequence of Kabat CDR-L1 of the humanized 14G8 antibody version 3

(Hu14G8VLv3).

SEQ ID NO: 81 sets forth the amino acid sequence of humanized 14G8 heavy chain version 1.

5 SEQ ID NO: 82 sets forth the amino acid sequence of humanized 14G8 heavy chain version 2.

SEQ ID NO: 83 sets forth the amino acid sequence of humanized 14G8 heavy chain version 3.

10 SEQ ID NO: 84 sets forth the amino acid sequence of humanized 14G8 light chain version 1.

SEQ ID NO: 85 sets forth the amino acid sequence of humanized 14G8 light chain version 2.

SEQ ID NO: 86 sets forth the amino acid sequence of humanized 14G8 light chain version 3.

15 SEQ ID NO: 87 sets forth the nucleic acid sequence of the heavy chain variable region of the mouse 14G8 antibody with signal peptide.

SEQ ID NO: 88 sets forth the amino acid sequence of the heavy chain variable region of the mouse 14G8 antibody with signal peptide.

20 SEQ ID NO: 89 sets forth the nucleic acid sequence of the light chain variable region of the mouse 14G8 antibody with signal peptide.

25 SEQ ID NO: 90 sets forth the amino acid sequence of the light chain variable region of the mouse 14G8 antibody with signal peptide.

SEQ ID NO: 91 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 14G8 antibody version 1 (Hu14G8VHv1).

30 SEQ ID NO: 92 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 14G8 antibody version 2 (Hu14G8VHv2).

SEQ ID NO: 93 sets forth the nucleic acid sequence of the heavy chain variable region of the humanized 14G8 antibody version 3 (Hu14G8VHv3).

35 SEQ ID NO: 94 sets forth the nucleic acid sequence of the light chain variable region of the humanized 14G8 antibody version 1 (Hu14G8VLv1).

40 SEQ ID NO: 95 sets forth the nucleic acid sequence of the light chain variable region of the humanized 14G8 antibody version 2 (Hu14G8VLv2).

SEQ ID NO: 96 sets forth the nucleic acid sequence of the light chain variable region of the humanized 14G8 antibody version 3 (Hu14G8VLv3).

45 SEQ ID NO: 97 sets forth the amino acid sequence of the mouse 14G8 heavy chain variable region signal peptide.

SEQ ID NO: 98 sets forth the nucleic acid sequence of the mouse 14G8 heavy chain variable region signal peptide.

50 SEQ ID NO: 99 sets forth the amino acid sequence of the mouse 14G8 light chain variable region signal peptide.

SEQ ID NO: 100 sets forth the nucleic acid sequence of the mouse 14G8 light chain variable region signal peptide.

SEQ ID NO: 101 sets forth the amino acid sequence of an exemplary human IgG1 heavy chain constant region.

55 SEQ ID NO: 102 sets forth the amino acid sequence of an exemplary human IgG1 heavy chain constant region of the IgG1 G1m3 allotype.

SEQ ID NO: 103 sets forth the amino acid sequence of an exemplary human IgG1 heavy chain constant region of

the IgG1 G1m3 allotype.

SEQ ID NO: 104 sets forth the amino acid sequence of an exemplary human kappa light chain constant region having an N-terminal arginine.

SEQ ID NO: 105 sets forth the amino acid sequence of an exemplary human kappa light chain constant region without an N-terminal arginine.

SEQ ID NO: 106 sets forth the nucleic acid sequence of an exemplary heavy chain constant region of the G1m3 allotype.

SEQ ID NO: 107 sets forth the nucleic acid sequence of an exemplary light chain constant region having an N-terminal arginine.

SEQ ID NO: 108 sets forth the nucleic acid sequence of an exemplary light chain constant region without an N-terminal arginine.

SEQ ID NO: 109 sets forth the amino acid sequence of human transthyretin set forth in accession number P02766.1 (UniProt).

SEQ ID NO: 110 sets forth the amino acid sequence of human transthyretin set forth in accession number AAB35639.1 (GenBank).

SEQ ID NO: 111 sets forth the amino acid sequence of human transthyretin set forth in accession number AAB35640.1 (GenBank).

SEQ ID NO: 112 sets forth the amino acid sequence of human transthyretin set forth in accession number and ABI63351.1 (GenBank).

SEQ ID NO: 113 sets forth the amino acid sequence of residues 89-97 of human transthyretin.

SEQ ID NO: 114 sets forth the amino acid sequence of a potential transthyretin immunogen.

SEQ ID NO: 115 sets forth the amino acid sequence of a potential transthyretin immunogen.

SEQ ID NO: 116 sets forth the amino acid sequence of a potential transthyretin immunogen.

SEQ ID NO: 117 sets forth the amino acid sequence of composite Chothia-Kabat CDR-H1 of the mouse 9D5 antibody.

SEQ ID NO: 118 sets forth the amino acid sequence of composite Chothia-Kabat CDR-H1 of the mouse 14G8 antibody.

## DEFINITIONS

**[0043]** Monoclonal antibodies or other biological entities are typically provided in isolated form. This means that an antibody or other biological entity is typically at least 50% w/w pure of interfering proteins and other contaminants arising from its production or purification but does not exclude the possibility that the monoclonal antibody is combined with an excess of pharmaceutically acceptable carrier(s) or other vehicle intended to facilitate its use. Sometimes monoclonal antibodies are at least 60%, 70%, 80%, 90%, 95% or 99% w/w pure of interfering proteins and contaminants from production or purification. Often an isolated monoclonal antibody or other biological entity is the predominant macromolecular species remaining after its purification.

**[0044]** Specific binding of an antibody to its target antigen means an affinity of at least  $10^6$ ,  $10^7$ ,  $10^8$ ,  $10^9$ , or  $10^{10}$  M<sup>-1</sup>. Specific binding is detectably higher in magnitude and distinguishable from non-specific binding occurring to at least one unrelated target. Specific binding can be the result of formation of bonds between particular functional groups or particular spatial fit (e.g., lock and key type) whereas nonspecific binding is usually the result of van der Waals forces. Specific binding does not however necessarily imply that an antibody binds one and only one target.

**[0045]** The basic antibody structural unit is a tetramer of subunits. Each tetramer includes two identical pairs of polypeptide chains, each pair having one "light" (about 25 kDa) and one "heavy" chain (about 50-70 kDa). The amino-terminal

portion of each chain includes a variable region of about 100 to 110 or more amino acids primarily responsible for antigen recognition. This variable region is initially expressed linked to a cleavable signal peptide. The variable region without the signal peptide is sometimes referred to as a mature variable region. Thus, for example, a light chain mature variable region means a light chain variable region without the light chain signal peptide. The carboxy-terminal portion of each chain defines a constant region primarily responsible for effector function.

**[0046]** Light chains are classified as either kappa or lambda. Heavy chains are classified as gamma, mu, alpha, delta, or epsilon, and define the antibody's isotype as IgG, IgM, IgA, IgD and IgE, respectively. Within light and heavy chains, the variable and constant regions are joined by a "J" region of about 12 or more amino acids, with the heavy chain also including a "D" region of about 10 or more amino acids. See *generally*, Fundamental Immunology, Paul, W., ed., 2nd ed. Raven Press, N.Y., 1989, Ch. 7.

**[0047]** An immunoglobulin light or heavy chain variable region (also referred to herein as a "light chain variable domain" ("VL domain") or "heavy chain variable domain" ("VH domain"), respectively) consists of a "framework" region interrupted by three "complementarity determining regions" or "CDRs." The framework regions serve to align the CDRs for specific binding to an epitope of an antigen. The CDRs include the amino acid residues of an antibody that are primarily responsible for antigen binding. From amino-terminus to carboxyl-terminus, both VL and VH domains comprise the following framework (FR) and CDR regions: FR1, CDR1, FR2, CDR2, FR3, CDR3, and FR4. CDRs 1, 2, and 3 of a VL domain are also referred to herein, respectively, as CDR-L1, CDR-L2, and CDR-L3; CDRs 1, 2, and 3 of a VH domain are also referred to herein, respectively, as CDR-H1, CDR-H2, and CDR-H3.

**[0048]** The assignment of amino acids to each VL and VH domain is in accordance with any conventional definition of CDRs. Conventional definitions include, the Kabat definition (Kabat, Sequences of Proteins of Immunological Interest (National Institutes of Health, Bethesda, MD, 1987 and 1991), The Chothia definition (Chothia & Lesk, J. Mol. Biol. 196:901-917, 1987; Chothia et al., Nature 342:878-883, 1989); a composite of Chothia Kabat CDR in which CDR-H1 is a composite of Chothia and Kabat CDRs; the AbM definition used by Oxford Molecular's antibody modelling software; and, the contact definition of Martin et al (bioinfo.org.uk/abs) (see Table 1). Kabat provides a widely used numbering convention (Kabat numbering) in which corresponding residues between different heavy chains or between different light chains are assigned the same number. When an antibody is said to comprise CDRs by a certain definition of CDRs (e.g., Kabat) that definition specifies the minimum number of CDR residues present in the antibody (*i.e.*, the Kabat CDRs). It does not exclude that other residues falling within another conventional CDR definition but outside the specified definition are also present. For example, an antibody comprising CDRs defined by Kabat includes among other possibilities, an antibody in which the CDRs contain Kabat CDR residues and no other CDR residues, and an antibody in which CDR H1 is a composite Chothia-Kabat CDR H1 and other CDRs contain Kabat CDR residues and no additional CDR residues based on other definitions.

Table 1

Conventional Definitions of CDRs Using Kabat Numbering					
Loop	Kabat	Chothia	Composite of Chothia & Kabat	AbM	Contact
L1	L24--L34	L24--L34	L24--L34	L24--L34	L30--L36
L2	L50--L56	L50--L56	L50--L56	L50--L56	L46--L55
L3	L89--L97	L89--L97	L89--L97	L89--L97	L89--L96
H1	H31--H35B	H26--H32..H34*	H26--H35B*	H26--H35B	H30--H35B
H2	H50--H65	H52--H56	H50--H65	H50--H58	H47--H58
H3	H95--H102	H95--H102	H95--H102	H95--H102	H93--H101

\*CDR-H1 by Chothia can end at H32, H33, or H34 (depending on the length of the loop). This is because the Kabat numbering scheme places insertions of extra residues at 35A and 35B, whereas Chothia numbering places them at 31A and 31B. If neither H35A nor H35B (Kabat numbering) is present, the Chothia CDR-H1 loop ends at H32. If only H35A is present, it ends at H33. If both H35A and H35B are present, it ends at H34.

**[0049]** The term "antibody" includes intact antibodies and binding fragments thereof. Typically, fragments compete with the intact antibody from which they were derived for specific binding to the target including separate heavy chains, light chains Fab, Fab', F(ab')<sub>2</sub>, F(ab)c, Dabs, nanobodies, and Fv. Fragments can be produced by recombinant DNA techniques, or by enzymatic or chemical separation of intact immunoglobulins. The term "antibody" also includes a bispecific antibody and/or a humanized antibody. A bispecific or bifunctional antibody is an artificial hybrid antibody having two different heavy/light chain pairs and two different binding sites (see, e.g., Songsivilai and Lachmann, Clin.

Exp. Immunol., 79:315-321 (1990); Kostelny et al., J. Immunol., 148:1547-53 (1992)). In some bispecific antibodies, the two different heavy/light chain pairs include a humanized 9D5 heavy chain/light chain pair and a heavy chain/light chain pair specific for a different epitope on transthyretin than that bound by 9D5. In some bispecific antibodies, the two different heavy/light chain pairs include a humanized 14G8 heavy chain/light chain pair and a heavy chain/light chain pair specific for a different epitope on transthyretin than that bound by 14G8.

**[0050]** In some bispecific antibodies, one heavy chain/light chain pair is a humanized 9D5 or 14G8 antibody as further disclosed below and the other heavy chain/light chain pair is from an antibody that binds to a receptor expressed on the blood brain barrier, such as an insulin receptor, an insulin-like growth factor (IGF) receptor, a leptin receptor, or a lipoprotein receptor, or a transferrin receptor (Friden et al., Proc. Natl. Acad. Sci. USA 88:4771-4775, 1991; Friden et al., Science 259:373-377, 1993). Such a bispecific antibody can be transferred cross the blood brain barrier by receptor-mediated transcytosis. Brain uptake of the bispecific antibody can be further enhanced by engineering the bi-specific antibody to reduce its affinity to the blood brain barrier receptor. Reduced affinity for the receptor resulted in a broader distribution in the brain (see, e.g., Atwal et al., Sci. Trans. Med. 3, 84ra43, 2011; Yu et al., Sci. Trans. Med. 3, 84ra44, 2011).

**[0051]** Exemplary bispecific antibodies can also be: (1) a dual-variable-domain antibody (DVD-Ig), where each light chain and heavy chain contains two variable domains in tandem through a short peptide linkage (Wu et al., Generation and Characterization of a Dual Variable Domain Immunoglobulin (DVD-Ig™) Molecule, In: Antibody Engineering, Springer Berlin Heidelberg (2010)); (2) a Tandab, which is a fusion of two single chain diabodies resulting in a tetravalent bispecific antibody that has two binding sites for each of the target antigens; (3) a flexibody, which is a combination of scFvs with a diabody resulting in a multivalent molecule; (4) a so-called "dock and lock" molecule, based on the "dimerization and docking domain" in Protein Kinase A, which, when applied to Fabs, can yield a trivalent bispecific binding protein consisting of two identical Fab fragments linked to a different Fab fragment; or (5) a so-called Scorpion molecule, comprising, e.g., two scFvs fused to both termini of a human Fc-region. Examples of platforms useful for preparing bispecific antibodies include BiTE (Micromet), DART (MacroGenics), Fcab and Mab2 (F-star), Fc-engineered IgG1 (Xencor) or DuoBody (based on Fab arm exchange, Genmab).

**[0052]** The term "epitope" refers to a site on an antigen to which an antibody binds. An epitope can be formed from contiguous amino acids or noncontiguous amino acids juxtaposed by tertiary folding of one or more proteins. Epitopes formed from contiguous amino acids (also known as linear epitopes) are typically retained on exposure to denaturing solvents whereas epitopes formed by tertiary folding (also known as conformational epitopes) are typically lost on treatment with denaturing solvents. An epitope typically includes at least 3, and more usually, at least 5 or 8-10 amino acids in a unique spatial conformation. Methods of determining spatial conformation of epitopes include, for example, x-ray crystallography and 2-dimensional nuclear magnetic resonance. See, e.g., Epitope Mapping Protocols, in Methods in Molecular Biology, Vol. 66, Glenn E. Morris, Ed. (1996). The epitope can be linear, such as an epitope of, for example, 2-5, 3-5, 3-9, or 5-9 contiguous amino acids from SEQ ID NO: 109. The epitope can also be a conformational epitope including, for example, two or more non-contiguous segments of amino acids within residues 89-97 of SEQ ID NO: 109.

**[0053]** If an antibody is said to bind to an epitope within amino acids 89-97 of transthyretin (TTR), for example, what is meant is that the epitope is within the recited range of amino acids including those defining the outer-limits of the range. It does not necessarily mean that every amino acid within the range constitutes part of the epitope. Thus, for example, an epitope within amino acids 89-97 of TTR may consist of amino acids 89-97, 89-96, 90-96, 91-96, 92-96, 93-96, 94-96, 89-96, 89-95, 89-94, 89-93, 89-92 or 89-93, among other linear segments of SEQ ID NO: 113, or in the case of conformational epitopes, non-contiguous segments of amino acids of SEQ ID NO: 113.

**[0054]** Antibodies that recognize the same or overlapping epitopes can be identified in a simple immunoassay showing the ability of one antibody to compete with the binding of another antibody to a target antigen. The epitope of an antibody can also be defined X-ray crystallography of the antibody bound to its antigen to identify contact residues. Alternatively, two antibodies have the same epitope if all amino acid mutations in the antigen that reduce or eliminate binding of one antibody reduce or eliminate binding of the other. Two antibodies have overlapping epitopes if some amino acid mutations that reduce or eliminate binding of one antibody reduce or eliminate binding of the other.

**[0055]** Competition between antibodies is determined by an assay in which an antibody under test inhibits specific binding of a reference antibody to a common antigen (see, e.g., Junghans *et al.*, *Cancer Res.* 50:1495, 1990). A test antibody competes with a reference antibody if an excess of a test antibody (e.g., at least 2x, 5x, 10x, 20x or 100x) inhibits binding of the reference antibody by at least 50% as measured in a competitive binding assay. Some test antibodies inhibit binding of the reference antibody by at least 75%, 90% or 99%. Antibodies identified by competition assay (competing antibodies) include antibodies binding to the same epitope as the reference antibody and antibodies binding to an adjacent epitope sufficiently proximal to the epitope bound by the reference antibody for steric hindrance to occur.

**[0056]** The term "native" with respect to the structure transthyretin (TTR) refers to the normal folded structure of TTR in its properly functioning state (*i.e.*, a TTR tetramer). As TTR is a tetramer in its natively folded form, non-native forms of TTR include, for example, misfolded TTR tetramers, TTR monomers, aggregated forms of TTR, and fibril forms of

TTR. Non-native forms of TTR can include molecules comprising wild-type TTR amino acid sequences or mutations.

**[0057]** The term "misfolded" with respect to TTR refers to the secondary and tertiary structure of a TTR polypeptide monomer or multimer, and indicates that the polypeptide has adopted a conformation that is not normal for that protein in its properly functioning state. Although TTR misfolding can be caused by mutations in the protein (e.g., deletion, substitution, or addition), wild-type TTR proteins can also be misfolded in diseases, exposing specific epitopes.

**[0058]** The term "pharmaceutically acceptable" means that the carrier, diluent, excipient, or auxiliary is compatible with the other ingredients of the formulation and not substantially deleterious to the recipient thereof.

**[0059]** The term "patient" includes human and other mammalian subjects that receive either prophylactic or therapeutic treatment.

**[0060]** An individual is at increased risk of a disease if the subject has at least one known risk-factor (e.g., genetic, biochemical, family history, and situational exposure) placing individuals with that risk factor at a statistically significant greater risk of developing the disease than individuals without the risk factor.

**[0061]** The term "biological sample" refers to a sample of biological material within or obtainable from a biological source, for example a human or mammalian subject. Such samples can be organs, organelles, tissues, sections of tissues, bodily fluids, peripheral blood, blood plasma, blood serum, cells, molecules such as proteins and peptides, and any parts or combinations derived therefrom. The term biological sample can also encompass any material derived by processing the sample. Derived material can include cells or their progeny. Processing of the biological sample may involve one or more of filtration, distillation, extraction, concentration, fixation, inactivation of interfering components, and the like.

**[0062]** The term "control sample" refers to a biological sample not known or suspected to include monomeric, misfolded, aggregated, or fibril forms of transthyretin (TTR), such as in TTR amyloid deposits. Control samples can be obtained from individuals not afflicted with a TTR amyloidosis or a specifically chosen type of TTR amyloidosis. Alternatively, control samples can be obtained from patients afflicted with TTR amyloidosis or a specifically chosen type of TTR amyloidosis. Such samples can be obtained at the same time as a biological sample thought to comprise the TTR amyloidosis or on a different occasion. A biological sample and a control sample can both be obtained from the same tissue (e.g., a tissue section containing both TTR amyloid deposits and surrounding normal tissue). Preferably, control samples consist essentially or entirely of tissue free of TTR amyloid deposits and can be used in comparison to a biological sample thought to comprise TTR amyloid deposits. Preferably, the tissue in the control sample is the same type as the tissue in the biological sample (e.g., cardiomyocytes in the heart).

**[0063]** The term "disease" refers to any abnormal condition that impairs physiological function. The term is used broadly to encompass any disorder, illness, abnormality, pathology, sickness, condition, or syndrome in which physiological function is impaired, irrespective of the nature of the etiology.

**[0064]** The term "symptom" refers to a subjective evidence of a disease, such as altered gait, as perceived by the subject. A "sign" refers to objective evidence of a disease as observed by a physician.

**[0065]** For purposes of classifying amino acids substitutions as conservative or nonconservative, amino acids are grouped as follows: Group I (hydrophobic side chains): met, ala, val, leu, ile; Group II (neutral hydrophilic side chains): cys, ser, thr; Group III (acidic side chains): asp, glu; Group IV (basic side chains): asn, gln, his, lys, arg; Group V (residues influencing chain orientation): gly, pro; and Group VI (aromatic side chains): trp, tyr, phe. Conservative substitutions involve substitutions between amino acids in the same class. Nonconservative substitutions constitute exchanging a member of one of these classes for a member of another.

**[0066]** Percentage sequence identities are determined with antibody sequences maximally aligned by the Kabat numbering convention. After alignment, if a subject antibody region (e.g., the entire mature variable region of a heavy or light chain) is being compared with the same region of a reference antibody, the percentage sequence identity between the subject and reference antibody regions is the number of positions occupied by the same amino acid in both the subject and reference antibody region divided by the total number of aligned positions of the two regions, with gaps not counted, multiplied by 100 to convert to percentage.

**[0067]** Compositions or methods "comprising" or "including" one or more recited elements may include other elements not specifically recited. For example, a composition that "comprises" or "includes" an antibody may contain the antibody alone or in combination with other ingredients.

**[0068]** Designation of a range of values includes all integers within or defining the range, and all subranges defined by integers within the range.

**[0069]** Unless otherwise apparent from the context, the term "about" encompasses values within a standard margin of error of measurement (e.g., SEM) of a stated value.

**[0070]** Statistical significance means  $p \leq 0.05$ .

**[0071]** The singular forms of the articles "a," "an," and "the" include plural references unless the context clearly dictates otherwise. For example, the term "a compound" or "at least one compound" can include a plurality of compounds, including mixtures thereof.

**DETAILED DESCRIPTION****I. General**

5 **[0072]** The disclosure provides antibodies that specifically bind to residues 89-97 of transthyretin (TTR). The antibodies have the capacity to bind to monomeric, misfolded, aggregated, or fibril forms of TTR. The antibodies can be used for treating or effecting prophylaxis of diseases or disorders associated with TTR accumulation or accumulation of TTR deposits (e.g., TTR amyloidosis). The antibodies can also be used for diagnosing TTR amyloidosis and inhibiting or reducing aggregation of TTR, among other applications.

**II. Target Molecules**

15 **[0073]** Transthyretin (TTR) is a 127-amino acid, 55 kDa serum and cerebrospinal fluid transport protein primarily synthesized by the liver. It has also been referred to as prealbumin, thyroxine binding prealbumin, ATTR, and TBPA. In its native state, TTR exists as a tetramer. In homozygotes, the tetramers comprise identical 127-amino-acid beta-sheet-rich subunits. In heterozygotes, the TTR tetramers are made up of variant and/or wild-type subunits, typically combined in a statistical fashion.

20 **[0074]** The established function of TTR in the blood is to transport holo-retinol binding protein. Although TTR is the major carrier of thyroxine ( $T_4$ ) in the blood of rodents, utilizing binding sites that are orthogonal to those used for holo-retinol binding protein, the  $T_4$  binding sites are effectively unoccupied in humans.

25 **[0075]** TTR is one of at least thirty different human proteins whose extracellular misfolding and/or misassembly (amyloidogenesis) into a spectrum of aggregate structures is thought to cause degenerative diseases referred to as amyloid diseases. TTR undergoes conformational changes in order to become amyloidogenic. Partial unfolding exposes stretches of largely uncharged hydrophobic residues in an extended conformation that efficiently misassemble into largely unstructured spherical aggregates that ultimately undergo conformation conversion into cross-beta sheet amyloid structures.

30 **[0076]** Unless otherwise apparent from context, reference to transthyretin (TTR) or its fragments or domains includes the natural human amino acid sequences including isoforms, mutants, and allelic variants thereof. Exemplary TTR polypeptide sequences are designated by Accession Numbers P02766.1 (UniProt), AAB35639.1 (GenBank), AAB35640.1 (GenBank), and AB163351.1 (GenBank) (SEQ ID NOS: 109-112, respectively). Residues are numbered according to Swiss Prot P02766.1, with the first amino acid of the mature protein (*i.e.*, not including the 20 amino acid signal sequence) designated residue 1. In any other TTR protein, residues are numbered according to the corresponding residues in P02766.1 on maximum alignment.

**III. Transthyretin Amyloidosis**

35 **[0077]** Transthyretin (TTR) amyloidosis is a systemic disorder characterized by pathogenic, misfolded TTR and the extracellular deposition of amyloid fibrils composed of TTR. TTR amyloidosis is generally caused by destabilization of the native TTR tetramer form (due to environmental or genetic conditions), leading to dissociation, misfolding, and aggregation of TTR into amyloid fibrils that accumulate in various organs and tissues, causing progressive dysfunction. See, e.g., Almeida and Saraiva, FEBS Letters 586:2891-2896 (2012); Ando et al., Orphanet Journal of Rare Diseases 8:31 (2013).

40 **[0078]** In humans, both wild-type TTR tetramers and mixed tetramers comprised of mutant and wild-type subunits can dissociate, misfold, and aggregate, with the process of amyloidogenesis leading to the degeneration of post-mitotic tissue. Thus, TTR amyloidoses encompass diseases caused by pathogenic misfolded TTR resulting from mutations in TTR or resulting from non-mutated, misfolded TTR.

45 **[0079]** For example, senile systemic amyloidosis (SSA) and senile cardiac amyloidosis (SCA) are age-related types of amyloidosis that result from the deposition of wild-type TTR amyloid outside and within the cardiomyocytes of the heart. TTR amyloidosis is also the most common form of hereditary (familial) amyloidosis, which is caused by mutations that destabilize the TTR protein. The TTR amyloidoses associated with point mutations in the TTR gene include familial amyloid polyneuropathy (FAP), familial amyloid cardiomyopathy (FAC), and the rare central nervous system selective amyloidosis (CNSA). Patients with hereditary (familial) TTR amyloidosis are almost always heterozygotes, meaning that the TTR tetramers are composed of mutant and/or wild-type TTR subunits, generally statistically distributed. Hereditary (familial) versions of TTR amyloidosis are generally autosomal dominant and are typically earlier onset than the sporadic diseases (SSA and SCA).

50 **[0080]** There are over 100 mutations in the gene encoding TTR that have been implicated in the autosomal dominant disorders FAP and FAC. See, e.g., US 2014/0056904; Saraiva, Hum. Mutat. 17(6):493-503 (2001); Damas and Saraiva, J. Struct. Biol. 130:290-299; Dwulet and Benson, Biochem. Biophys. Res. Commun. 114:657-662 (1983). These amyloid-

causing mutations are distributed throughout the entire molecule of TTR. Generally, the more destabilizing the mutant subunits are to the TTR tetramer structure, the earlier the onset of amyloid disease. The pathogenic potential of a TTR variant is generally determined by a combination of its instability and its cellular secretion efficiency. The initial pathology caused by some TTR variants comes from their selective destruction of cardiac tissue, whereas that from other TTR variants comes from compromising the peripheral and autonomic nervous system. The tissue damage caused by TTR amyloidogenesis appear to stem largely from the toxicity of small, diffusible TTR aggregates, although accumulation of extracellular amyloid may contribute and almost certainly compromises organ structure in the late stages of the TTR amyloidosis.

**[0081]** TTR amyloidosis presents in many different forms, with considerable phenotypic variation across individuals and geographic locations. For example, TTR amyloidosis can present as a progressive, axonal sensory autonomic and motor neuropathy. TTR amyloidosis can also present as an infiltrative cardiomyopathy.

**[0082]** The age at onset of disease-related symptoms varies between the second and ninth decades of life, with great variations across different populations. The multisystem involvement of TTR amyloidosis is a clue to its diagnosis. For example, TTR amyloidosis diagnosis is considered when one or several of the following are present: (1) family history of neuropathic disease, especially associated with heart failure; (2) neuropathic pain or progressive sensory disturbances of unknown etiology; (3) carpal tunnel syndrome without obvious cause, particularly if it is bilateral and requires surgical release; (4) gastrointestinal motility disturbances or autonomic nerve dysfunction of unknown etiology (e.g., erectile dysfunction, orthostatic hypotension, neurogenic bladder); (5) cardiac disease characterized by thickened ventricular walls in the absence of hypertension; (6) advanced atrio-ventricular block of unknown origin, particularly when accompanied by a thickened heart; and (6) vitreous body inclusions of the cotton-wool type. See Ando et al., *Orphanet Journal of Rare Diseases* 8:31 (2013). Other symptoms can include, for example, polyneuropathy, sensory loss, pain, weakness in lower limbs, dyshidrosis, diarrhea, constipation, weight loss, and urinary incontinence/retention.

**[0083]** Diagnosis of TTR amyloidosis typically relies on target organ biopsies, followed by histological staining of the excised tissue with the amyloid-specific dye, Congo red. If a positive test for amyloid is observed, immunohistochemical staining for TTR is subsequently performed to ensure that the precursor protein responsible for amyloid formation is indeed TTR. For familial forms of the diseases, demonstration of a mutation in the gene encoding TTR is then needed before diagnosis can be made. This can be accomplished, for example, through isoelectric focusing electrophoresis, polymerase chain reaction, or laser dissection/liquid chromatography-tandem mass spectrometry. See, e.g., US 2014/0056904; Ruberg and Berk, *Circulation* 126:1286-1300 (2012); Ando et al., *Orphanet Journal of Rare Diseases* 8:31 (2013).

#### IV. Antibodies

##### A. Binding Specificity and Functional Properties

**[0084]** The disclosure provides monoclonal antibodies binding to transthyretin (TTR) protein, more specifically, to epitopes within amino acid residues 89-97 (SEQ ID NO: 113) of TTR. Such epitopes are buried in the native TTR tetramer and exposed in monomeric, misfolded, aggregated, or fibril forms of TTR.

**[0085]** Antibodies designated 9D5 and 14G8 are two such exemplary mouse antibodies. Unless otherwise apparent from context, reference to 9D5 or 14G8 should be understood as referring to any of the mouse, chimeric, veneered, and humanized forms of these antibodies. These antibodies specifically bind within amino acid residues 89-97 (SEQ ID NO: 113) of TTR. These antibodies are further characterized by their ability to bind to monomeric, misfolded, aggregated, or fibril forms of TTR but not to native tetrameric forms of TTR. In addition, these antibodies are characterized by their immunoreactivity on TTR-mediated amyloidosis cardiac tissue but not on healthy cardiac tissue. Ability to bind to specific proteins or fragments thereof may be demonstrated using exemplary assay formats provided in the examples.

**[0086]** Some antibodies bind to the same or overlapping epitope as an antibody designated 9D5 or 14G8. The sequences of the heavy and light chain mature variable regions of these antibodies are designated SEQ ID NOS: 1 and 16 (9D5), and 61 and 70 (14G8), respectively. Other antibodies having such a binding specificity can be produced by immunizing mice with TTR, or a portion thereof including the desired epitope (e.g., SEQ ID NO: 113), and screening resulting antibodies for binding to monomeric TTR or a peptide comprising SEQ ID NO: 113, optionally in competition with an antibody having the variable regions of mouse 9D5 or 14G8 (IgG1 kappa). Fragments of TTR including the desired epitope can be linked to a carrier that helps elicit an antibody response to the fragment and/or be combined with an adjuvant that helps elicit such a response. Such antibodies can be screened for differential binding to wild-type, monomeric versions of TTR or a fragment thereof (e.g., SEQ ID NO: 113) compared with mutants of specified residues. Screening against such mutants more precisely defines the binding specificity to allow identification of antibodies whose binding is inhibited by mutagenesis of particular residues and which are likely to share the functional properties of other exemplified antibodies. The mutations can be systematic replacement substitution with alanine (or serine if an alanine is present already) one residue at a time, or more broadly spaced intervals, throughout the target or throughout a section

thereof in which an epitope is known to reside. If the same set of mutations significantly reduces the binding of two antibodies, the two antibodies bind the same epitope.

**[0087]** Antibodies having the binding specificity of a selected murine antibody (e.g., 9D5 or 14G8) can also be produced using a variant of the phage display method. See Winter, WO 92/20791. This method is particularly suitable for producing human antibodies. In this method, either the heavy or light chain variable region of the selected murine antibody is used as a starting material. If, for example, a light chain variable region is selected as the starting material, a phage library is constructed in which members display the same light chain variable region (*i.e.*, the murine starting material) and a different heavy chain variable region. The heavy chain variable regions can for example be obtained from a library of rearranged human heavy chain variable regions. A phage showing strong specific binding (e.g., at least  $10^8$  and preferably at least  $10^9$  M<sup>-1</sup>) for monomeric TTR or a fragment thereof (e.g., amino acid residues 89-97) is selected. The heavy chain variable region from this phage then serves as a starting material for constructing a further phage library. In this library, each phage displays the same heavy chain variable region (*i.e.*, the region identified from the first display library) and a different light chain variable region. The light chain variable regions can be obtained for example from a library of rearranged human variable light chain regions. Again, phage showing strong specific binding for monomeric TTR or a fragment thereof (e.g., amino acid residues 89-97) are selected. The resulting antibodies usually have the same or similar epitope specificity as the murine starting material. Kabat CDRs of the heavy chain of 9D5 are designated SEQ ID NOS: 13-15, respectively, and Kabat CDRs of the light chain of 9D5 are designated SEQ ID NOS: 24-26, respectively. A composite Chothia-Kabat CDR-H1 of 9D5 is designated SEQ ID NO: 117. Kabat CDRs of the heavy chain of 14G8 are designated SEQ ID NOS: 67-69, respectively, and Kabat CDRs of the light chain of 14G8 are designated SEQ ID NOS: 77-79, respectively. A composite Chothia-Kabat CDR-H1 of 14G8 is designated SEQ ID NO: 118. A variant of CDR-L1 of 14G8 is designated SEQ ID NO: 80.

**[0088]** Other antibodies can be obtained by mutagenesis of cDNA encoding the heavy and light chains of an exemplary antibody, such as 9D5 or 14G8. Monoclonal antibodies that are at least 70%, 80%, 90%, 95%, 96%, 97%, 98%, or 99% identical to 9D5 or 14G8 in amino acid sequence of the mature heavy and/or light chain variable regions and maintain its functional properties, and/or which differ from the respective antibody by a small number of functionally inconsequential amino acid substitutions (e.g., conservative substitutions), deletions, or insertions are also included in the invention. Monoclonal antibodies having at least one or all six CDR(s) as defined by any conventional definition, but preferably Kabat, that are 90%, 95%, 99% or 100% identical to corresponding CDRs of 9D5 or 14G8 are also included.

**[0089]** The disclosure also provides antibodies having some or all (*e.g.*, 3, 4, 5, and 6) CDRs entirely or substantially from 9D5 or 14G8. Such antibodies can include a heavy chain variable region that has at least two, and usually all three, CDRs entirely or substantially from the heavy chain variable region of 9D5 or 14G8 and/or a light chain variable region having at least two, and usually all three, CDRs entirely or substantially from the light chain variable region of 9D5 or 14G8. The antibodies can include both heavy and light chains. A CDR is substantially from a corresponding 9D5 or 14G8 CDR when it contains no more than 4, 3, 2, or 1 substitutions, insertions, or deletions, except that CDR-H2 (when defined by Kabat) can have no more than 6, 5, 4, 3, 2, or 1 substitutions, insertions, or deletions. Such antibodies can have at least 70%, 80%, 90%, 95%, 96%, 97%, 98%, or 99% identity to 9D5 or 14G8 in the amino acid sequence of the mature heavy and/or light chain variable regions and maintain their functional properties, and/or differ from 9D5 or 14G8 by a small number of functionally inconsequential amino acid substitutions (e.g., conservative substitutions), deletions, or insertions.

**[0090]** Some antibodies identified by such assays can bind to monomeric, misfolded, aggregated, or fibril forms of TTR but not to native tetrameric forms of TTR, as described in the examples or otherwise. Likewise, some antibodies are immunoreactive on TTR-mediated amyloidosis tissue but not on healthy tissue.

**[0091]** Some antibodies can inhibit or reduce aggregation of TTR, inhibit or reduce TTR fibril formation, reduce or clear TTR deposits or aggregated TTR, or stabilize non-toxic conformations of TTR in an animal model or clinical trial. Some antibodies can treat, effect prophylaxis of, or delay the onset of a TTR amyloidosis as shown in an animal model or clinical trial. Exemplary animal models for testing activity against a TTR amyloidosis include those described in Kohno et al., *Am. J. Path.* 150(4): 1497-1508 (1997); Teng et al., *Laboratory Investigations* 81:385-396 (2001); Wakasugi et al., *Proc. Japan Acad.* 63B:344-347 (1987); Shimada et al., *Mol. Biol. Med.* 6:333-343 (1989); Nagata et al., *J. Biochem.* 117:169-175 (1995); Sousa et al., *Am. J. Path.* 161:1935-1948 (2002); and Santos et al., *Neurobiology of Aging* 31:280-289 (2010).

## **B. Non-Human Antibodies**

**[0092]** The production of other non-human antibodies, *e.g.*, murine, guinea pig, primate, rabbit or rat, against monomeric TTR or a fragment thereof (e.g., amino acid residues 89-97) can be accomplished by, for example, immunizing the animal with TTR or a fragment thereof. See Harlow & Lane, *Antibodies, A Laboratory Manual* (CSHP NY, 1988). Such an immunogen can be obtained from a natural source, by peptide synthesis, or by recombinant expression. Optionally, the immunogen can be administered fused or otherwise complexed with a carrier protein. Optionally, the immunogen

can be administered with an adjuvant. Several types of adjuvant can be used as described below. Complete Freund's adjuvant followed by incomplete adjuvant is preferred for immunization of laboratory animals. Rabbits or guinea pigs are typically used for making polyclonal antibodies. Mice are typically used for making monoclonal antibodies. Antibodies are screened for specific binding to monomeric TTR or an epitope within TTR (e.g., an epitope comprising one or more of amino acid residues 89-97). Such screening can be accomplished by determining binding of an antibody to a collection of monomeric TTR variants, such as TTR variants containing amino acid residues 89-97 or mutations within these residues, and determining which TTR variants bind to the antibody. Binding can be assessed, for example, by Western blot, FACS or ELISA.

### C. Humanized Antibodies

**[0093]** A humanized antibody is a genetically engineered antibody in which CDRs from a non-human "donor" antibody are grafted into human "acceptor" antibody sequences (see, e.g., Queen, US 5,530,101 and 5,585,089; Winter, US 5,225,539; Carter, US 6,407,213; Adair, US 5,859,205; and Foote, US 6,881,557). The acceptor antibody sequences can be, for example, a mature human antibody sequence, a composite of such sequences, a consensus sequence of human antibody sequences, or a germline region sequence. Thus, a humanized antibody is an antibody having at least three, four, five or all CDRs entirely or substantially from a donor antibody and variable region framework sequences and constant regions, if present, entirely or substantially from human antibody sequences. Similarly a humanized heavy chain has at least one, two and usually all three CDRs entirely or substantially from a donor antibody heavy chain, and a heavy chain variable region framework sequence and heavy chain constant region, if present, substantially from human heavy chain variable region framework and constant region sequences. Similarly a humanized light chain has at least one, two and usually all three CDRs entirely or substantially from a donor antibody light chain, and a light chain variable region framework sequence and light chain constant region, if present, substantially from human light chain variable region framework and constant region sequences. Other than nanobodies and dAbs, a humanized antibody comprises a humanized heavy chain and a humanized light chain. A CDR in a humanized antibody is substantially from a corresponding CDR in a non-human antibody when at least 85%, 90%, 95% or 100% of corresponding residues (as defined by any conventional definition but preferably defined by Kabat) are identical between the respective CDRs. The variable region framework sequences of an antibody chain or the constant region of an antibody chain are substantially from a human variable region framework sequence or human constant region respectively when at least 85%, 90%, 95% or 100% of corresponding residues defined by any conventional definition but preferably defined by Kabat are identical.

**[0094]** Although humanized antibodies often incorporate all six CDRs (defined by any conventional definition but preferably as defined by Kabat) from a mouse antibody, they can also be made with less than all CDRs (e.g., at least 3, 4, or 5 CDRs) from a mouse antibody (e.g., Pascalis et al., J. Immunol. 169:3076, 2002; Vajdos et al., J. of Mol. Biol., 320:415-428, 2002; Iwahashi et al., Mol. Immunol. 36:1079-1091, 1999; Tamura et al., J. Immunol., 164:1432-1441, 2000).

**[0095]** In some antibodies only part of the CDRs, namely the subset of CDR residues required for binding, termed the SDRs, are needed to retain binding in a humanized antibody. CDR residues not contacting antigen and not in the SDRs can be identified based on previous studies (for example residues H60-H65 in CDR H2 are often not required), from regions of Kabat CDRs lying outside Chothia hypervariable loops (Chothia, J. Mol. Biol. 196:901, 1987), by molecular modeling and/or empirically, or as described in Gonzales et al., Mol. Immunol. 41: 863, 2004. In such humanized antibodies at positions in which one or more donor CDR residues is absent or in which an entire donor CDR is omitted, the amino acid occupying the position can be an amino acid occupying the corresponding position (by Kabat numbering) in the acceptor antibody sequence. The number of such substitutions of acceptor for donor amino acids in the CDRs to include reflects a balance of competing considerations. Such substitutions are potentially advantageous in decreasing the number of mouse amino acids in a humanized antibody and consequently decreasing potential immunogenicity. However, substitutions can also cause changes of affinity, and significant reductions in affinity are preferably avoided. Positions for substitution within CDRs and amino acids to substitute can also be selected empirically.

**[0096]** The human acceptor antibody sequences can optionally be selected from among the many known human antibody sequences to provide a high degree of sequence identity (e.g., 65-85% identity) between a human acceptor sequence variable region frameworks and corresponding variable region frameworks of a donor antibody chain.

**[0097]** Examples of acceptor sequences for the heavy chain are the human mature heavy chain variable regions with NCBI accession codes BAC02114 and AAX82494.1 (SEQ ID NOS: 3 and 4) and heavy chain variable regions of human Kabat subgroup 3. BAC02114 shares the same canonical form as mouse 9D5 heavy chain. Other examples of acceptor sequences for the heavy chain are the human mature heavy chain variable regions with NCBI accession codes AAD30410.1 and AAX82494.1 (SEQ ID NOS: 63 and 4, respectively) and heavy chain variable regions of human Kabat subgroup 1. AAD30410.1 and AAX82494.1 include two CDRs having the same canonical form as mouse 14G8 heavy chain. Examples of acceptor sequences for the light chain are the human mature light chain variable region with NCBI accession code ABC66952 (SEQ ID NO: 18) and light chain variable regions of human Kabat subgroup 3. ABC66952 includes two CDRs having the same canonical form as mouse 9D5 light chain. Other examples of acceptor sequences

for the light chain are the human mature light chain variable regions with NCBI accession codes ABA71374.1 and ABC66952.1 (SEQ ID NOS: 72 and 73, respectively) and light chain variable regions of human Kabat subgroup 2. ABA71374.1 and ABC66952.1 have the same canonical form as mouse 14G8 light chain.

**[0098]** If more than one human acceptor antibody sequence is selected, a composite or hybrid of those acceptors can be used, and the amino acids used at different positions in the humanized light chain and heavy chain variable regions can be taken from any of the human acceptor antibody sequences used. For example, the human mature heavy chain variable regions with NCBI accession codes BAC02114 and AAX82494.1 (SEQ ID NOS: 3 and 4) were used as acceptor sequences for humanization of the 9D5 mature heavy chain variable region. Examples of positions in which these two acceptors differ include positions H19 (R or K), H40 (A or T), H44 (G or R), H49 (S or A), H77 (S or T), H82a (N or S), H83 (R or K), H84 (A or S), and H89 (V or M). Humanized versions of the 9D5 heavy chain variable region can include either amino acid at any of these positions. Similarly, the human mature heavy chain variable regions with NCBI accession codes AAD30410.1 and AAX82494.1 (SEQ ID NOS: 63 and 4, respectively) were used as acceptor sequences for humanization of the 14G8 mature heavy chain variable region. Examples of positions in which these two acceptors differ include positions H82a (N or S), H83 (R or K), H84 (A or S), and H89 (V or M). Humanized versions of the 14G8 heavy chain variable region can include either amino acid at any of these positions. Similarly, the human mature light chain variable regions with NCBI accession codes ABA71374.1 and ABC66952.1 (SEQ ID NOS: 72 and 73, respectively) were used as acceptor sequences for humanization of the 14G8 mature light chain variable region. An example of a position in which these two acceptors differ is position L18 (S or P). Humanized versions of the 14G8 light chain variable region can include either amino acid at this position.

**[0099]** Certain amino acids from the human variable region framework residues can be selected for substitution based on their possible influence on CDR conformation and/or binding to antigen. Investigation of such possible influences is by modeling, examination of the characteristics of the amino acids at particular locations, or empirical observation of the effects of substitution or mutagenesis of particular amino acids.

**[0100]** For example, when an amino acid differs between a murine variable region framework residue and a selected human variable region framework residue, the human framework amino acid can be substituted by the equivalent framework amino acid from the mouse antibody when it is reasonably expected that the amino acid:

- (1) noncovalently binds antigen directly;
- (2) is adjacent to a CDR region or within a CDR as defined by Chothia but not Kabat;
- (3) otherwise interacts with a CDR region (*e.g.*, is within about 6 Å of a CDR region), (*e.g.*, identified by modeling the light or heavy chain on the solved structure of a homologous known immunoglobulin chain); or
- (4) is a residue participating in the VL-VH interface.

**[0101]** Framework residues from classes (1) through (3) as defined by Queen, US 5,530,101, are sometimes alternately referred to as canonical and vernier residues. Framework residues that help define the conformation of a CDR loop are sometimes referred to as canonical residues (Chothia & Lesk, *J. Mol. Biol.* 196:901-917 (1987); Thornton & Martin, *J. Mol. Biol.* 263:800-815 (1996)). Framework residues that support antigen-binding loop conformations and play a role in fine-tuning the fit of an antibody to antigen are sometimes referred to as vernier residues (Foote & Winter, *J. Mol. Biol.* 224:487-499 (1992)).

**[0102]** Other framework residues that are candidates for substitution are residues creating a potential glycosylation site. Still other candidates for substitution are acceptor human framework amino acids that are unusual for a human immunoglobulin at that position. These amino acids can be substituted with amino acids from the equivalent position of the mouse donor antibody or from the equivalent positions of more typical human immunoglobulins.

**[0103]** Exemplary humanized antibodies are humanized forms of the mouse 9D5 or 14G8 antibodies, designated Hu9D5 or Hu14G8, respectively. The mouse 9D5 antibody comprises mature heavy and light chain variable regions having amino acid sequences comprising SEQ ID NO: 1 and SEQ ID NO: 16, respectively. The disclosure provides eight exemplified humanized mature heavy chain variable regions: Hu9D5VHv1, Hu9D5VHv2, Hu9D5VHv2b, Hu9D5VHv3, Hu9D5VHv3b, Hu9D5VHv4, Hu9D5VHv4b, and Hu9D5VHv5 (SEQ ID NOS: 5-12, respectively). The disclosure further provides five exemplified human mature light chain variable regions: Hu9D5VLv1, Hu9D5VLv2, Hu9D5VLv3, Hu9D5VLv4, and Hu9D5VLv5 (SEQ ID NOS: 19-23, respectively). Figure 1 shows alignments of 9D5, mouse model antibodies, human acceptor antibodies, and various humanized antibodies.

**[0104]** The mouse 14G8 antibody comprises mature heavy and light chain variable regions having amino acid sequences comprising SEQ ID NO: 61 and SEQ ID NO: 70, respectively. The disclosure provides three exemplified humanized mature heavy chain variable regions: Hu14G8VHv1, Hu14G8VHv2, and Hu14G8VHv3 (SEQ ID NOS: 64-66, respectively). The disclosure further provides three exemplified human mature light chain variable regions: Hu14G8VLv1, Hu14G8VLv2, and Hu14G8VLv3 (SEQ ID NOS: 74-76, respectively). Figure 2 shows alignments of 14G8, mouse model antibodies, human acceptor antibodies, and various humanized antibodies.

**[0105]** For reasons such as possible influence on CDR conformation and/or binding to antigen, mediating interaction

between heavy and light chains, interaction with the constant region, being a site for desired or undesired post-translational modification, being an unusual residue for its position in a human variable region sequence and therefore potentially immunogenic, getting aggregation potential, and other reasons, the following 15 variable region framework positions were considered as candidates for substitutions in the eight exemplified Hu9D5 mature heavy chain variable regions and the five exemplified Hu9D5 mature light chain variable regions, as further specified in the examples: H42 (G42E), H47 (W47L), H69 (I69F), H82 (M82S), H82b (S82(b)L), H108 (T108L), L8 (P8A), L9 (L9P), L18 (P18S), L19 (A19V), L36 (Y36F), L39 (K39R), L60 (D60S), L70 (D70A), and L74 (K74R). Likewise, the following 11 variable region framework positions were considered as candidates for substitutions in the three exemplified Hu14G8 mature heavy chain variable regions and the three exemplified Hu14G8 mature light chain variable regions, as further specified in the examples: H1 (Q1E), H3 (Q3K), H47 (W47L), H105 (Q105T), L8 (P8A), L9 (L9P), L19 (A19V), L26 (N26S), L36 (Y36F), L60 (D60S), and L70 (D70A).

**[0106]** Here, as elsewhere, the first-mentioned residue is the residue of a humanized antibody formed by grafting Kabat CDRs or a composite Chothia-Kabat CDR in the case of CDR-H1 into a human acceptor framework (e.g., a composite or hybrid human acceptor framework), and the second-mentioned residue is a residue being considered for replacing such residue. Thus, within variable region frameworks, the first mentioned residue is human, and within CDRs, the first mentioned residue is mouse.

**[0107]** Exemplified Hu9D5 antibodies include any permutations or combinations of the exemplified mature heavy and light chain variable regions (e.g., VHv1/VLv1 or H1L1, VHv1/VLv2 or H1L2, VHv1/VLv3 or H1L3, VHv1/VLv4 or H1L4, VHv1/VLv5 or H1L5, VHv2/VLvl or H2L1, VHv2/VLv2 or H2L2, VHv2/VLv3 or H2L3, VHv2/VLv4 or H2L4, VHv2/VLv5 or H2L5, VHv2b/VLv1 or H2bL1, VHv2b/VLv2 or H2bL2, VHv2b/VLv3 or H2bL3, VHv2b/VLv4 or H2bL4, VHv2b/VLv5 or H2bL5, VHv3/VLvl or H3L1, VHv3/VLv2 or H3L2, VHv3/VLv3 or H3L3, VHv3/VLv4 or H3L4, VHv3/VLv5 or H3L5, VHv3b/VLv1 or H3bL1, VHv3b/VLv2 or H3bL2, VHv3b/VLv3 or H3bL3, VHv3b/VLv4 or H3bL4, VHv3b/VLv5 or H3bL5, VHv4/VLvl or H4L1, VHv4/VLv2 or H4L2, VHv4/VLv3 or H4L3, VHv4/VLv4 or H4L4, VHv4/VLv5 or H4L5, VHv4b/VLv1 or H4bL1, VHv4b/VLv2 or H4bL2, VHv4b/VLv3 or H4bL3, VHv4b/VLv4 or H4bL4, VHv4b/VLv5 or H4bL5, VHv5/VLv1 or H5L1, VHv5/VLv2 or H5L2, VHv5/VLv3 or H5L3, VHv5/VLv4 or H5L4, and VHv5/VLv5 or H5L5).

**[0108]** The disclosure provides variants of humanized 9D5 antibodies in which the humanized mature heavy chain variable region shows at least 90%, 95%, 96%, 97%, 98%, or 99% identity to a humanized Hu9D5VHv4b (SEQ ID NO: 11) and the humanized mature light chain variable region shows at least 90%, 95%, 96%, 97%, 98%, or 99% identity to a Hu9D5VLv1 (SEQ ID NO: 19). In some such antibodies, at least 1, 2, or all 3 of the backmutations or other mutations in Hu9D5 H4bL1 are retained. The disclosure also provides variants of the other exemplified humanized 9D5 antibodies. Such variants have mature light and heavy chain variable regions showing at least 90%, 95%, 96%, 97%, 98%, or 99% sequence identity to the mature light and heavy chain variable regions of the exemplified humanized 9D5 H1L1, H1L2, H1L3, H1L4, H1L5, H2L1, H2L2, H2L3, H2L4, H2L5, H2bL1, H2bL2, H2bL3, H2bL4, H2bL5, H3L1, H3L2, H3L3, H3L4, H3L5, H3bL1, H3bL2, H3bL3, H3bL4, H3bL5, H4L1, H4L2, H4L3, H4L4, H4L5, H4bL1, H4bL2, H4bL3, H4bL4, H4bL5, H5L1, H5L2, H5L3, H5L4, or H5L5 antibodies.

**[0109]** In some antibodies, at least one of positions H42, H47, H69, H82, H82b, and H108 in the Vh region is occupied by E, L, F, S, L, and L, respectively. In some antibodies, positions H47, H69, and H82 in the Vh region are occupied by L, F, and S, respectively, as in Hu9D5VHv1. In some antibodies, positions H47, H69, H82, and H82b in the Vh region are occupied by L, F, S, and L, respectively, as in Hu9D5VHv2. In some antibodies, positions H42, H47, and H108 in the Vh region are occupied by E, L, and L, respectively, as in Hu9D5VHv2b. In some antibodies, positions H69, H82, and H82b in the Vh region are occupied by F, S, and L, respectively, as in Hu9D5VHv3. In some antibodies, positions H47 and H108 in the Vh region are each occupied by L, as in Hu9D5VHv3b and Hu9D5Vh4b. In some antibodies, positions H82 and H82b in the Vh region are occupied by S and L, respectively, as in Hu9D5VHv4. In some antibodies, positions H42, H47, and H82b in the Vh region are occupied by E, L, and L, respectively, as in Hu9D5VHv5. In some antibodies, at least one of positions L8, L9, L18, L19, L36, L39, L60, L70, and L74 in the Vk region is occupied by A, P, S, V, F, R, S, A, and R, respectively. In some antibodies, position L36 in the Vk region is occupied by F, as in Hu9D5VLvl. In some antibodies, position L60 in the Vk region is occupied by S, as in Hu9D5VLv3. In some antibodies, positions L8, L9, L19, L36, L39, L60, L70, and L74 in the Vk region are occupied by A, P, V, F, R, S, A, and R, respectively, as in Hu9D5VLv4. In some antibodies, positions L8, L9, L18, L19, L36, L39, L60, L70, and L74 in the Vk region are occupied by A, P, S, V, F, R, S, A, and R, respectively, as in Hu9D5VLv5. The CDR regions of such humanized antibodies can be identical or substantially identical to the CDR regions of the 9D5 mouse donor antibody or any of the above exemplified humanized 9D5 antibodies. The CDR regions can be defined by any conventional definition (e.g., Chothia, or composite of Chothia and Kabat) but are preferably as defined by Kabat.

**[0110]** Variable regions framework positions are in accordance with Kabat numbering unless otherwise stated. Other such variants typically differ from the sequences of the exemplified Hu9D5 heavy and light chains by a small number (e.g., typically no more than 1, 2, 3, 5, 10, or 15) of replacements, deletions or insertions. Such differences are usually in the framework but can also occur in the CDRs.

**[0111]** Exemplified Hu14G8 antibodies include any permutations or combinations of the exemplified mature heavy

and light chain variable regions (e.g., VHv1/VLv1 or H1L1, VHv1/VLv2 or H1L2, VHv1/VLv3 or H1L3, VHv2/VLv1 or H2L1, VHv2/VLv2 or H2L2, VHv2/VLv3 or H2L3, VHv3/VLv1 or H3L1, VHv3/VLv2 or H3L2, and VHv3/VLv3 or H3L3).

[0112] The disclosure provides variants of humanized 14G8 antibodies in which the humanized mature heavy chain variable region shows at least 90%, 95%, 96%, 97%, 98%, or 99% identity to Hu14G8VHv2 (Hu14G8 H2) (SEQ ID NO: 65) and the humanized mature light chain variable region shows at least 90%, 95%, 96%, 97%, 98%, or 99% identity to Hu14G8VLv3 (Hu14G8 L3) (SEQ ID NO: 76). In some such antibodies, at least 1, 2, 3, 4, or all 5 of the backmutations or other mutations in Hu14G8H2L3 are retained. The disclosure also provides variants of the other exemplified humanized 14G8 antibodies. Such variants have mature light and heavy chain variable regions showing at least 90%, 95%, 96%, 97%, 98%, or 99% sequence identity to the mature light and heavy chain variable regions of the exemplified humanized 14G8 H1L1, H1L2, H1L3, H2L1, H2L2, H2L3, H3L1, H3L2, or H3L3 antibodies.

[0113] In some antibodies, at least one of positions H1 and H47 in the Vh region is occupied by E and L, respectively. In some antibodies, positions H1 and H47 in the Vh region are occupied by E and L, respectively, as in Hu14G8VHv2 and Hu14G8VHv3. In some antibodies, at least one of positions H3 and H105 in the Vh region is occupied by K and T, respectively. In some antibodies, positions H3 and H105 in the Vh region are occupied by K and T, respectively, as in Hu14G8VHv1. In some antibodies, position L36 in the Vk region is occupied by F, as in Hu14G8VLv2. In some antibodies, at least one of positions L8, L9, L19, L26, L60, and L70 in the Vk region is occupied by A, P, V, S, S, and A, respectively. In some antibodies, positions L8, L9, L19, and L70 in the Vk region are occupied by A, P, V, and A, respectively, as in Hu14G8VLv1. In some antibodies, positions L26 and L60 in the Vk region are each occupied by S, as in Hu14G8VLv3. The CDR regions of such humanized antibodies can be identical or substantially identical to the CDR regions of the 14G8 mouse donor antibody. The CDR regions can be defined by any conventional definition (e.g., Chothia, or composite of Chothia and Kabat) but are preferably as defined by Kabat.

[0114] Variable regions framework positions are in accordance with Kabat numbering unless otherwise stated. Other such variants typically differ from the sequences of the exemplified Hu14G8 heavy and light chains by a small number (e.g., typically no more than 1, 2, 3, 5, 10, or 15) of replacements, deletions or insertions. Such differences are usually in the framework but can also occur in the CDRs

[0115] A possibility for additional variation in humanized 9D5 or 14G8 variants is additional backmutations in the variable region frameworks. Many of the framework residues not in contact with the CDRs in the humanized mAb can accommodate substitutions of amino acids from the corresponding positions of the donor mouse mAb or other mouse or human antibodies, and even many potential CDR-contact residues are also amenable to substitution. Even amino acids within the CDRs may be altered, for example, with residues found at the corresponding position of the human acceptor sequence used to supply variable region frameworks. In addition, alternate human acceptor sequences can be used, for example, for the heavy and/or light chain. If different acceptor sequences are used, one or more of the backmutations recommended above may not be performed because the corresponding donor and acceptor residues are already the same without backmutations.

[0116] Preferably, replacements or backmutations in humanized 9D5 or 14G8 variants (whether or not conservative) have no substantial effect on the binding affinity or potency of the humanized mAb, that is, its ability to bind to monomeric TTR (e.g., the potency in some or all of the assays described in the present examples of the variant humanized 9D5 or 14G8 antibody is essentially the same, *i.e.*, within experimental error, as that of murine 9D5 or 14G8 antibody).

#### **D. Chimeric and Veneered Antibodies**

[0117] The disclosure further provides chimeric and veneered forms of non-human antibodies, particularly the 9D5 or 14G8 antibodies of the examples.

[0118] A chimeric antibody is an antibody in which the mature variable regions of light and heavy chains of a non-human antibody (e.g., a mouse) are combined with human light and heavy chain constant regions. Such antibodies substantially or entirely retain the binding specificity of the mouse antibody, and are about two-thirds human sequence.

[0119] A veneered antibody is a type of humanized antibody that retains some and usually all of the CDRs and some of the non-human variable region framework residues of a non-human antibody but replaces other variable region framework residues that may contribute to B- or T-cell epitopes, for example exposed residues (Padlan, *Mol. Immunol.* 28:489, 1991) with residues from the corresponding positions of a human antibody sequence. The result is an antibody in which the CDRs are entirely or substantially from a non-human antibody and the variable region frameworks of the non-human antibody are made more human-like by the substitutions. Veneered forms of the 9D5 or 14G8 antibody are included in the disclosure.

#### **E. Human Antibodies**

[0120] Human antibodies against monomeric TTR or a fragment thereof (e.g., amino acid residues 89-97 (SEQ ID NO: 113) of TTR) are provided by a variety of techniques described below. Some human antibodies are selected by

competitive binding experiments, by the phage display method of Winter, above, or otherwise, to have the same epitope specificity as a particular mouse antibody, such as one of the mouse monoclonal antibodies described in the examples. Human antibodies can also be screened for particular epitope specificity by using only a fragment of TTR, such as a TTR variant containing only amino acid residues 89-97 of TTR, as the target antigen, and/or by screening antibodies against a collection of TTR variants, such as TTR variants containing various mutations within amino acid residues 89-97 of TTR.

**[0121]** Methods for producing human antibodies include the trioma method of Oestberg et al., Hybridoma 2:361-367 (1983); Oestberg, U.S. Patent No. 4,634,664; and Engleman et al., US Patent 4,634,666, use of transgenic mice including human immunoglobulin genes (see, e.g., Lonberg et al., WO93/12227 (1993); US 5,877,397; US 5,874,299; US 5,814,318; US 5,789,650; US 5,770,429; US 5,661,016; US 5,633,425; US 5,625,126; US 5,569,825; US 5,545,806; Neuberger, Nat. Biotechnol. 14:826 (1996); and Kucherlapati, WO 91/10741 (1991)) and phage display methods (see, e.g., Dower et al., WO 91/17271; McCafferty et al., WO 92/01047; US 5,877,218; US 5,871,907; US 5,858,657; US 5,837,242; US 5,733,743; and US 5,565,332).

#### **F. Selection of Constant Region**

**[0122]** The heavy and light chain variable regions of chimeric, veneered or humanized antibodies can be linked to at least a portion of a human constant region. The choice of constant region depends, in part, whether antibody-dependent cell-mediated cytotoxicity, antibody dependent cellular phagocytosis and/or complement dependent cytotoxicity are desired. For example, human isotopes IgG1 and IgG3 have complement-dependent cytotoxicity and human isotopes IgG2 and IgG4 do not. Human IgG1 and IgG3 also induce stronger cell mediated effector functions than human IgG2 and IgG4. Light chain constant regions can be lambda or kappa.

**[0123]** One or several amino acids at the amino or carboxy terminus of the light and/or heavy chain, such as the C-terminal lysine of the heavy chain, may be missing or derivatized in a proportion or all of the molecules. Substitutions can be made in the constant regions to reduce or increase effector function such as complement-mediated cytotoxicity or ADCC (see, e.g., Winter et al., US Patent No. 5,624,821; Tso et al., US Patent No. 5,834,597; and Lazar et al., Proc. Natl. Acad. Sci. USA 103:4005, 2006), or to prolong half-life in humans (see, e.g., Hinton et al., J. Biol. Chem. 279:6213, 2004). Exemplary substitutions include a Gln at position 250 and/or a Leu at position 428 (EU numbering is used in this paragraph for the constant region) for increasing the half-life of an antibody. Substitution at any or all of positions 234, 235, 236 and/or 237 reduce affinity for Fc $\gamma$  receptors, particularly Fc $\gamma$ RI receptor (see, e.g., US 6,624,821). An alanine substitution at positions 234, 235, and 237 of human IgG1 can be used for reducing effector functions. Some antibodies have alanine substitution at positions 234, 235 and 237 of human IgG1 for reducing effector functions. Optionally, positions 234, 236 and/or 237 in human IgG2 are substituted with alanine and position 235 with glutamine (see, e.g., US 5,624,821). In some antibodies, a mutation at one or more of positions 241, 264, 265, 270, 296, 297, 322, 329, and 331 by EU numbering of human IgG1 is used. In some antibodies, a mutation at one or more of positions 318, 320, and 322 by EU numbering of human IgG1 is used. In some antibodies, positions 234 and/or 235 are substituted with alanine and/or position 329 is substituted with glycine. In some antibodies, positions 234 and 235 are substituted with alanine, such as in SEQ ID NO: 102. In some antibodies, the isotype is human IgG2 or IgG4.

**[0124]** An exemplary human light chain kappa constant region has the amino acid sequence of SEQ ID NO: 104. The N-terminal arginine of SEQ ID NO: 104 can be omitted, in which case light chain kappa constant region has the amino acid sequence of SEQ ID NO: 105. An exemplary human IgG1 heavy chain constant region has the amino acid sequence of SEQ ID NO: 101 (with or without the C-terminal lysine). Antibodies can be expressed as tetramers containing two light and two heavy chains, as separate heavy chains, light chains, as Fab, Fab', F(ab')<sub>2</sub>, and Fv, or as single chain antibodies in which heavy and light chain mature variable domains are linked through a spacer.

**[0125]** Human constant regions show allotypic variation and isoallotypic variation between different individuals, that is, the constant regions can differ in different individuals at one or more polymorphic positions. Isoallotypes differ from allotypes in that sera recognizing an isoallotype bind to a non-polymorphic region of a one or more other isotopes. Thus, for example, another heavy chain constant region is of IgG1 G1m3 allotype and has the amino acid sequence of SEQ ID NO: 103. Another heavy chain constant region of the IgG1 G1m3 allotype has the amino acid sequence of SEQ ID NO: 102 (with or without the C-terminal lysine). Reference to a human constant region includes a constant region with any natural allotype or any permutation of residues occupying positions in natural allotypes.

#### **G. Expression of Recombinant Antibodies**

**[0126]** A number of methods are known for producing chimeric and humanized antibodies using an antibody-expressing cell line (e.g., hybridoma). For example, the immunoglobulin variable regions of antibodies can be cloned and sequenced using well known methods. In one method, the heavy chain variable VH region is cloned by RT-PCR using mRNA prepared from hybridoma cells. Consensus primers are employed to the VH region leader peptide encompassing the

translation initiation codon as the 5' primer and a g2b constant regions specific 3' primer. Exemplary primers are described in U.S. patent publication US 2005/0009150 by Schenk *et al.* (hereinafter "Schenk"). The sequences from multiple, independently derived clones can be compared to ensure no changes are introduced during amplification. The sequence of the VH region can also be determined or confirmed by sequencing a VH fragment obtained by 5' RACE RT-PCR methodology and the 3' g2b specific primer.

**[0127]** The light chain variable VL region can be cloned in an analogous manner. In one approach, a consensus primer set is designed for amplification of VL regions using a 5' primer designed to hybridize to the VL region encompassing the translation initiation codon and a 3' primer specific for the Ck region downstream of the V-J joining region. In a second approach, 5'RACE RT-PCR methodology is employed to clone a VL encoding cDNA. Exemplary primers are described in Schenk, *supra*. The cloned sequences are then combined with sequences encoding human (or other non-human species) constant regions. Exemplary sequences encoding human constant regions include SEQ ID NO: 106, which encodes a human IgG1 constant region (SEQ ID NO: 103), and SEQ ID NOS: 107 and 108, which encode human kappa light chain constant regions (SEQ ID NOS: 104 and 105, respectively).

**[0128]** In one approach, the heavy and light chain variable regions are re-engineered to encode splice donor sequences downstream of the respective VDJ or VJ junctions and are cloned into a mammalian expression vector, such as pCMV-hyl for the heavy chain and pCMV-Mcl for the light chain. These vectors encode human  $\gamma$ 1 and Ck constant regions as exonic fragments downstream of the inserted variable region cassette. Following sequence verification, the heavy chain and light chain expression vectors can be co-transfected into CHO cells to produce chimeric antibodies. Conditioned media is collected 48 hours post-transfection and assayed by western blot analysis for antibody production or ELISA for antigen binding. The chimeric antibodies are humanized as described above.

**[0129]** Chimeric, veneered, humanized, and human antibodies are typically produced by recombinant expression. Recombinant polynucleotide constructs typically include an expression control sequence operably linked to the coding sequences of antibody chains, including naturally associated or heterologous expression control elements, such as a promoter. The expression control sequences can be promoter systems in vectors capable of transforming or transfecting eukaryotic or prokaryotic host cells. Once the vector has been incorporated into the appropriate host, the host is maintained under conditions suitable for high level expression of the nucleotide sequences and the collection and purification of the crossreacting antibodies.

**[0130]** These expression vectors are typically replicable in the host organisms either as episomes or as an integral part of the host chromosomal DNA. Commonly, expression vectors contain selection markers, *e.g.*, ampicillin resistance or hygromycin resistance, to permit detection of those cells transformed with the desired DNA sequences.

**[0131]** *E. coli* is one prokaryotic host useful for expressing antibodies, particularly antibody fragments. Microbes, such as yeast, are also useful for expression. *Saccharomyces* is a yeast host with suitable vectors having expression control sequences, an origin of replication, termination sequences, and the like as desired. Typical promoters include 3-phosphoglycerate kinase and other glycolytic enzymes. Inducible yeast promoters include, among others, promoters from alcohol dehydrogenase, isocytochrome C, and enzymes responsible for maltose and galactose utilization.

**[0132]** Mammalian cells can be used for expressing nucleotide segments encoding immunoglobulins or fragments thereof. See Winnacker, *From Genes to Clones*, (VCH Publishers, NY, 1987). A number of suitable host cell lines capable of secreting intact heterologous proteins have been developed, and include CHO cell lines, various COS cell lines, HeLa cells, HEK293 cells, L cells, and non-antibody-producing myelomas including Sp2/0 and NS0. The cells can be nonhuman. Expression vectors for these cells can include expression control sequences, such as an origin of replication, a promoter, an enhancer (Queen *et al.*, *Immunol. Rev.* 89:49 (1986)), and necessary processing information sites, such as ribosome binding sites, RNA splice sites, polyadenylation sites, and transcriptional terminator sequences. Expression control sequences can include promoters derived from endogenous genes, cytomegalovirus, SV40, adenovirus, bovine papillomavirus, and the like. See Co *et al.*, *J. Immunol.* 148:1149 (1992).

**[0133]** Alternatively, antibody coding sequences can be incorporated in transgenes for introduction into the genome of a transgenic animal and subsequent expression in the milk of the transgenic animal (*see, e.g.*, U.S. Pat. No. 5,741,957; U.S. Pat. No. 5,304,489; and U.S. Pat. No. 5,849,992). Suitable transgenes include coding sequences for light and/or heavy chains operably linked with a promoter and enhancer from a mammary gland specific gene, such as casein or beta lactoglobulin.

**[0134]** The vectors containing the DNA segments of interest can be transferred into the host cell by methods depending on the type of cellular host. For example, calcium chloride transfection is commonly utilized for prokaryotic cells, whereas calcium phosphate treatment, electroporation, lipofection, biolistics, or viral-based transfection can be used for other cellular hosts. Other methods used to transform mammalian cells include the use of polybrene, protoplast fusion, liposomes, electroporation, and microinjection. For production of transgenic animals, transgenes can be microinjected into fertilized oocytes or can be incorporated into the genome of embryonic stem cells, and the nuclei of such cells transferred into enucleated oocytes.

**[0135]** Having introduced vector(s) encoding antibody heavy and light chains into cell culture, cell pools can be screened for growth productivity and product quality in serum-free media. Top-producing cell pools can then be subjected of FACS-

based single-cell cloning to generate monoclonal lines. Specific productivities above 50 pg or 100 pg per cell per day, which correspond to product titers of greater than 7.5 g/L culture, can be used. Antibodies produced by single cell clones can also be tested for turbidity, filtration properties, PAGE, IEF, UV scan, HP-SEC, carbohydrate-oligosaccharide mapping, mass spectrometry, and binding assay, such as ELISA or Biacore. A selected clone can then be banked in multiple vials and stored frozen for subsequent use.

**[0136]** Once expressed, antibodies can be purified according to standard procedures of the art, including protein A capture, HPLC purification, column chromatography, gel electrophoresis and the like (*see generally*, Scopes, Protein Purification (Springer-Verlag, NY, 1982)).

**[0137]** Methodology for commercial production of antibodies can be employed, including codon optimization, selection of promoters, selection of transcription elements, selection of terminators, serum-free single cell cloning, cell banking, use of selection markers for amplification of copy number, CHO terminator, or improvement of protein titers (*see, e.g.*, US 5,786,464; US 6,114,148; US 6,063,598; US 7,569,339; WO2004/050884; WO2008/012142; WO2008/012142; WO2005/019442; WO2008/107388; WO2009/027471; and US 5,888,809).

#### 15 **H. Antibody Screening Assays**

**[0138]** Antibodies can be subject to several screens including binding assays, functional screens, screens in animal models of diseases associated with TTR deposits, and clinical trials. Binding assays test for specific binding and, optionally, affinity and epitope specificity to monomeric TTR or a fragment thereof. For example, binding assays can screen for antibodies that bind to amino acid residues 89-97 (SEQ ID NO: 113) of TTR, which is an epitope that is buried in the native TTR tetramer and exposed in monomeric, misfolded, aggregated, or fibril forms of TTR. Antibodies can also be screened for the ability to bind pre-fibrillar, non-native conformations of TTR and TTR amyloid fibrils but not native TTR conformations. For example, antibodies can be screened for the ability to bind to monomeric forms of TTR created by dissociation or disaggregation of native tetrameric TTR, and can be counter-screened against native tetrameric TTR, as described in the examples or otherwise. Likewise, antibodies can also be screened for their immunoreactivity on TTR-mediated amyloidosis tissue but not on healthy tissue. Such screens are sometimes performed in competition with an exemplary antibody, such as an antibody having the variable regions of 9D5 or 14G8 (IgG1 kappa isotype). Optionally, either the antibody or TTR target is immobilized in such assay.

**[0139]** Functional assays can be performed in cellular models including cells naturally expressing TTR or transfected with DNA encoding TTR or a fragment thereof. Suitable cells include cells derived from cardiac tissue or other tissues affected by TTR amyloidogenesis. Cells can be screened for reduced levels of monomeric, misfolded, aggregated, or fibril forms of TTR (*e.g.*, by Western blotting or immunoprecipitation of cell extracts or supernatants) or reduced toxicity attributable to monomeric, misfolded, aggregated, or fibril forms of TTR. For example, antibodies can be tested for the ability to inhibit or reduce aggregation of TTR, inhibit or reduce TTR fibril formation, reduce TTR deposits, clear aggregated TTR, or stabilize non-toxic conformations of TTR.

**[0140]** Other functional assays can be performed in solution, such as testing whether an antibody is capable of disrupting or reducing TTR fibril formation when monomeric TTR or misfolded TTR intermediates in solution are contacted with the antibody. The extent of fibril formation can be probed by turbidity measurements, for example, at 400 nm on a UV-visible spectrometer equipped with a temperature control unit. Thioflavin-T can also be used to assess the extent of amyloid fibril formation. For example, a five-fold molar excess of Thioflavin-T can be added to TTR samples and left at room temperature for 30 minutes before measurements are taken. Thioflavin-T fluorescence can be monitored using a spectrofluorimeter. *See* US 2014/0056904.

**[0141]** Animal model screens test the ability of the antibody to therapeutically or prophylactically treat signs or symptoms in an animal model simulating a human disease associated with accumulation of TTR or TTR deposits. Such diseases include types of TTR amyloidosis, such as senile systemic amyloidosis (SSA), senile cardiac amyloidosis (SCA), familial amyloid polyneuropathy (FAP), familial amyloid cardiomyopathy (FAC), and central nervous system selective amyloidosis (CNSA). Suitable signs or symptoms that can be monitored include the presence and extent of amyloid deposits in various tissues, such as the gastrointestinal tract or heart. The extent of reduction of amyloid deposits can be determined by comparison with an appropriate control, such the level of TTR amyloid deposits in control animals that have received a control antibody (*e.g.*, an isotype matched control antibody), a placebo, or no treatment at all. An exemplary animal model for testing activity against a TTR amyloidosis is a mouse model carrying a null mutation at the endogenous mouse *Ttr* locus and the human mutant *TTR* gene comprising a V30M mutation that is associated with familial amyloidotic polyneuropathy. *See, e.g.*, Kohno et al., *Am. J. Path.* 150(4):1497-1508 (1997); Cardoso and Saraiva, *FASEB J* 20(2):234-239 (2006). Similar models also exist, including other models for familial versions of TTR amyloidosis and models for sporadic versions of TTR amyloidosis. *See, e.g.*, Teng et al., *Lab. Invest.* 81(3): 385-396 (2001); Ito and Maeda, *Mouse Models of Transthyretin Amyloidosis*, in *Recent Advances in Transthyretin Evolution, Structure, and Biological Functions*, pp. 261-280 (2009) (Springer Berlin Heidelberg). Transgenic animals can include a human TTR transgene, such as a TTR transgene with a mutation associated with TTR amyloidosis or a wild-type TTR transgene.

To facilitate testing in animal models, chimeric antibodies having a constant region appropriate for the animal model can be used (e.g., mouse-rat chimeras could be used for testing antibodies in rats). It can be concluded that a humanized version of an antibody will be effective if the corresponding mouse antibody or chimeric antibody is effective in an appropriate animal model and the humanized antibody has similar binding affinity (e.g., within experimental error, such as by a factor of 1.5, 2, or 3).

**[0142]** Clinical trials test for safety and efficacy in a human having a disease associated with TTR amyloidosis.

### ***I. Nucleic Acids***

**[0143]** The disclosure further provides nucleic acids encoding any of the heavy and light chains described above (e.g., SEQ ID NOS: 40, 42, 44-56, 87, 89, 91-96, and 106-108). Optionally, such nucleic acids further encode a signal peptide and can be expressed with the signal peptide linked to the constant region (e.g., signal peptides having amino acid sequences of SEQ ID NOS: 41 and 88 (heavy chain) and 43 and 90 (light chain) that can be encoded by SEQ ID NOS: 40 and 87, respectively (heavy chain) and 42 and 89, respectively (light chain)). Coding sequences of nucleic acids can be operably linked with regulatory sequences to ensure expression of the coding sequences, such as a promoter, enhancer, ribosome binding site, transcription termination signal, and the like. The nucleic acids encoding heavy and light chains can occur in isolated form or can be cloned into one or more vectors. The nucleic acids can be synthesized by, for example, solid state synthesis or PCR of overlapping oligonucleotides. Nucleic acids encoding heavy and light chains can be joined as one contiguous nucleic acid, e.g., within an expression vector, or can be separate, e.g., each cloned into its own expression vector.

### ***J. Conjugated Antibodies***

**[0144]** Conjugated antibodies that specifically bind to antigens exposed in pathogenic forms of TTR but not in native tetrameric forms of TTR, such as amino acid residues 89-97 (SEQ ID NO: 113) of TTR, are useful in detecting the presence of monomeric, misfolded, aggregated, or fibril forms of TTR; monitoring and evaluating the efficacy of therapeutic agents being used to treat patients diagnosed with a TTR amyloidosis; inhibiting or reducing aggregation of TTR; inhibiting or reducing TTR fibril formation; reducing or clearing TTR deposits; stabilizing non-toxic conformations of TTR; or treating or effecting prophylaxis of a TTR amyloidosis in a patient. For example, such antibodies can be conjugated with other therapeutic moieties, other proteins, other antibodies, and/or detectable labels. See WO 03/057838; US 8,455,622.

**[0145]** Conjugated therapeutic moieties can be any agent that can be used to treat, combat, ameliorate, prevent, or improve an unwanted condition or disease in a patient, such as a TTR amyloidosis. Therapeutic moieties can include, for example, immunomodulators or any biologically active agents that facilitate or enhance the activity of the antibody. An immunomodulator can be any agent that stimulates or inhibits the development or maintenance of an immunologic response. If such therapeutic moieties are coupled to an antibody specific for monomeric, misfolded, aggregated, or fibril forms of TTR, such as the antibodies described herein, the coupled therapeutic moieties will have a specific affinity for non-native, pathogenic forms of TTR over native tetrameric forms of TTR. Consequently, administration of the conjugated antibodies directly targets tissues comprising pathogenic forms of TTR with minimal damage to surrounding normal, healthy tissue. This can be particularly useful for therapeutic moieties that are too toxic to be administered on their own. In addition, smaller quantities of the therapeutic moieties can be used.

**[0146]** Examples of suitable therapeutic moieties include drugs that reduce levels of TTR, stabilize the native tetrameric structure of TTR, inhibit aggregation of TTR, disrupt TTR fibril or amyloid formation, or counteract cellular toxicity. See, e.g., Almeida and Saraiva, FEBS Letters 586:2891-2896 (2012); Saraiva, FEBS Letters 498:201-203 (2001); Ando et al., Orphanet Journal of Rare Diseases 8:31 (2013); Ruberg and Berk, Circulation 126:1286-1300 (2012); and Johnson et al., J. Mol. Biol. 421(2-3):185-203 (2012). For example, antibodies can be conjugated to tafamidis, diflunisal, ALN-TTR01, ALN-TTR02, ISIS-TTRRx, doxycycline (doxy), tauroursodeoxycholic acid (TUDCA), Doxy-TUDCA, epigallocatechin gallate (EGCG), curcumin, or resveratrol (3,5,4'-trihydroxystilbene). Other representative therapeutic moieties include other agents known to be useful for treatment, management, or amelioration of a TTR amyloidosis or symptoms of a TTR amyloidosis. See, e.g., Ando et al., Orphanet Journal of Rare Diseases 8:31 (2013) for common clinical symptoms of TTR amyloidosis and typical agents used to treat those symptoms.

**[0147]** Antibodies can also be coupled with other proteins. For example, antibodies can be coupled with Fynomers. Fynomers are small binding proteins (e.g., 7 kDa) derived from the human Fyn SH3 domain. They can be stable and soluble, and they can lack cysteine residues and disulfide bonds. Fynomers can be engineered to bind to target molecules with the same affinity and specificity as antibodies. They are suitable for creating multi-specific fusion proteins based on antibodies. For example, Fynomers can be fused to N-terminal and/or C-terminal ends of antibodies to create bi- and tri-specific FynomAbs with different architectures. Fynomers can be selected using Fynomer libraries through screening technologies using FACS, Biacore, and cell-based assays that allow efficient selection of Fynomers with optimal properties. Examples of Fynomers are disclosed in Grabulovski et al., J. Biol. Chem. 282:3196-3204 (2007); Bertschinger

et al., *Protein Eng. Des. Sel.* 20:57-68 (2007); Schlatter et al., *MAbs.* 4:497-508 (2011); Banner et al., *Acta. Crystallogr. D. Biol. Crystallogr.* 69(Pt6):1124-1137 (2013); and Brack et al., *Mol. Cancer Ther.* 13:2030-2039 (2014).

**[0148]** The antibodies disclosed herein can also be coupled or conjugated to one or more other antibodies (e.g., to form antibody heteroconjugates). Such other antibodies can bind to different epitopes within TTR or a portion thereof or can bind to a different target antigen.

**[0149]** Antibodies can also be coupled with a detectable label. Such antibodies can be used, for example, for diagnosing a TTR amyloidosis, for monitoring progression of a TTR amyloidosis, and/or for assessing efficacy of treatment. Such antibodies are particularly useful for performing such determinations in subjects having or being susceptible to a TTR amyloidosis, or in appropriate biological samples obtained from such subjects. Representative detectable labels that may be coupled or linked to an antibody include various enzymes, such as horseradish peroxidase, alkaline phosphatase, beta-galactosidase, or acetylcholinesterase; prosthetic groups, such as streptavidin/biotin and avidin/biotin; fluorescent materials, such as umbelliferone, fluorescein, fluorescein isothiocyanate, rhodamine, dichlorotriazinylamine fluorescein, dansyl chloride or phycoerythrin; luminescent materials, such as luminol; bioluminescent materials, such as luciferase, luciferin, and aequorin; radioactive materials, such as yttrium<sup>90</sup> (90Y), radiosilver-111, radiosilver-199, Bismuth<sup>213</sup>, iodine (<sup>131</sup>I, <sup>125</sup>I, <sup>123</sup>I, <sup>121</sup>I, ), carbon (<sup>14</sup>C), sulfur (<sup>35</sup>S), tritium (<sup>3</sup>H), indium (<sup>115</sup>In, <sup>113</sup>In, <sup>112</sup>In, <sup>111</sup>In, ), technetium (<sup>99</sup>Tc), thallium (<sup>201</sup>Tl), gallium (<sup>68</sup>Ga, <sup>67</sup>Ga), palladium (<sup>103</sup>Pd), molybdenum (<sup>99</sup>Mo), xenon (<sup>133</sup>Xe), fluorine (<sup>18</sup>F), <sup>153</sup>Sm, <sup>177</sup>Lu, <sup>159</sup>Gd, <sup>149</sup>Pm, <sup>140</sup>La, <sup>175</sup>Yb, <sup>166</sup>Ho, <sup>90</sup>Y, <sup>47</sup>Sc, <sup>186</sup>Re, <sup>188</sup>Re, <sup>142</sup>Pr, <sup>105</sup>Rh, <sup>97</sup>Ru, <sup>68</sup>Ge, <sup>57</sup>Co, <sup>65</sup>Zn, <sup>85</sup>Sr, <sup>32</sup>P, <sup>153</sup>Gd, <sup>169</sup>Yb, <sup>51</sup>Cr, <sup>54</sup>Mn, <sup>75</sup>Se, <sup>113</sup>Sn, and <sup>117</sup>Tm; positron emitting metals using various positron emission tomographies; nonradioactive paramagnetic metal ions; and molecules that are radiolabelled or conjugated to specific radioisotopes.

**[0150]** Linkage of radioisotopes to antibodies may be performed with conventional bifunctional chelates. For radiosilver-111 and radiosilver-199 linkage, sulfur-based linkers may be used. See Hazra et al., *Cell Biophys.* 24-25:1-7 (1994). Linkage of silver radioisotopes may involve reducing the immunoglobulin with ascorbic acid. For radioisotopes such as <sup>111</sup>In and <sup>90</sup>Y, ibritumomab tiuxetan can be used and will react with such isotopes to form <sup>111</sup>Inibritumomab tiuxetan and <sup>90</sup>Y-ibritumomab tiuxetan, respectively. See Witzig, *Cancer Chemother. Pharmacol.*, 48 Suppl 1:S91-S95 (2001).

**[0151]** Therapeutic moieties, other proteins, other antibodies, and/or detectable labels may be coupled or conjugated, directly or indirectly through an intermediate (e.g., a linker), to an antibody of the disclosure. See e.g., Aron et al., "Monoclonal Antibodies For Immunotargeting Of Drugs In Cancer Therapy," in *Monoclonal Antibodies And Cancer Therapy*, Reisfeld et al. (eds.), pp. 243-56 (Alan R. Liss, Inc. 1985); Hellstrom et al., "Antibodies For Drug Delivery," in *Controlled Drug Delivery (2nd Ed.)*, Robinson et al. (eds.), pp. 623-53 (Marcel Dekker, Inc. 1987); Thorpe, "Antibody Carriers Of Cytotoxic Agents In Cancer Therapy: A Review," in *Monoclonal Antibodies 84: Biological And Clinical Applications*, Pinchera et al. (eds.), pp. 475-506 (1985); "Analysis, Results, And Future Prospective Of The Therapeutic Use Of Radiolabeled Antibody In Cancer Therapy," in *Monoclonal Antibodies For Cancer Detection And Therapy*, Baldwin et al. (eds.), pp. 303-16 (Academic Press 1985); and Thorpe et al., *Immunol. Rev.*, 62:119-58 (1982). Suitable linkers include, for example, cleavable and non-cleavable linkers. Different linkers that release the coupled therapeutic moieties, proteins, antibodies, and/or detectable labels under acidic or reducing conditions, on exposure to specific proteases, or under other defined conditions can be employed.

## V. Therapeutic Applications

**[0152]** The above antibodies can be used for treating or effecting prophylaxis of a disease in a patient having or at risk for the disease mediated at least in part by transthyretin (TTR), and particularly by monomeric, misfolded, aggregated, or fibril forms of TTR. Although an understanding of mechanism is not required for practice, it is believed that any or all of the following mechanisms may contribute to treatment of TTR amyloidosis using the above antibodies: antibody-mediated inhibition of TTR aggregation and fibril formation, antibody-mediated stabilization of non-toxic conformations of TTR (e.g., tetrameric forms), or antibody-mediated clearance of aggregated TTR, oligomeric TTR, or monomeric TTR. Antibody-drug conjugates can have additional mechanisms of action determined by the conjugated moiety.

**[0153]** Antibodies are administered in an effective regime meaning a dosage, route of administration and frequency of administration that delays the onset, reduces the severity, inhibits further deterioration, and/or ameliorates at least one sign or symptom of a disorder being treated. If a patient is already suffering from a disorder, the regime can be referred to as a therapeutically effective regime. If the patient is at elevated risk of the disorder relative to the general population but is not yet experiencing symptoms, the regime can be referred to as a prophylactically effective regime. In some instances, therapeutic or prophylactic efficacy can be observed in an individual patient relative to historical controls or past experience in the same patient. In other instances, therapeutic or prophylactic efficacy can be demonstrated in a preclinical or clinical trial in a population of treated patients relative to a control population of untreated patients.

**[0154]** The frequency of administration depends on the half-life of the antibody in the circulation, the condition of the patient and the route of administration among other factors. The frequency can be daily, weekly, monthly, quarterly, or at irregular intervals in response to changes in the patient's condition or progression of the disorder being treated. An exemplary frequency for intravenous administration is between weekly and quarterly over a continuous course of treat-

ment, although more or less frequent dosing is also possible. For subcutaneous administration, an exemplary dosing frequency is daily to monthly, although more or less frequent dosing is also possible.

[0155] The number of dosages administered depends on whether the disorder is acute or chronic and the response of the disorder to the treatment. For acute disorders or acute exacerbations of a chronic disorder, between 1 and 10 doses are often sufficient. Sometimes a single bolus dose, optionally in divided form, is sufficient for an acute disorder or acute exacerbation of a chronic disorder. Treatment can be repeated for recurrence of an acute disorder or acute exacerbation. For chronic disorders, an antibody can be administered at regular intervals, e.g., weekly, fortnightly, monthly, quarterly, every six months for at least 1, 5 or 10 years, or the life of the patient.

## VI. Pharmaceutical Compositions and Methods of Use

[0156] Provided herein are several methods of diagnosing, monitoring, treating or effecting prophylaxis of diseases or conditions mediated at least in part by transthyretin (TTR), and particularly by monomeric, misfolded, aggregated, or fibril forms of TTR (e.g., TTR amyloidosis). Examples of such diseases include familial TTR amyloidoses, such as familial amyloid cardiomyopathy (FAC), familial amyloid polyneuropathy (FAP), or central nervous system selective amyloidosis (CNSA), and sporadic TTR amyloidoses, such as senile systemic amyloidosis (SSA) or senile cardiac amyloidosis (SCA). Antibodies described above can be incorporated into a pharmaceutical composition for use in such methods. In general, an antibody or pharmaceutical composition containing an antibody is administered to a subject in need thereof. Patients amenable to treatment include individuals at risk of TTR amyloidosis but not showing symptoms, as well as patients presently showing symptoms. Some patients can be treated during the prodromal stage of TTR amyloidosis.

[0157] The pharmaceutical compositions can be administered prophylactically to individuals who have a known genetic risk of TTR amyloidosis. Such individuals include those having relatives who have experienced such a disease, and those whose risk is determined by analysis of genetic or biochemical markers (e.g., mutations in TTR associated with TTR amyloidosis), including using the diagnostic methods provided herein. For example, there are over 100 mutations in the gene encoding TTR that have been implicated in TTR amyloidosis. See, e.g., US 2014/0056904; Saraiva, Hum. Mutat. 17(6):493-503 (2001); Damas and Saraiva, J. Struct. Biol. 130:290-299; Dwulet and Benson, Biochem. Biophys. Res. Commun. 114:657-662 (1983).

[0158] Individuals suffering from TTR amyloidosis can sometimes be recognized from the clinical manifestations of TTR amyloidosis, including one or more of the following: (1) family history of neuropathic disease, especially associated with heart failure; (2) neuropathic pain or progressive sensory disturbances of unknown etiology; (3) carpal tunnel syndrome without obvious cause, particularly if it is bilateral and requires surgical release; (4) gastrointestinal motility disturbances or autonomic nerve dysfunction of unknown etiology (e.g., erectile dysfunction, orthostatic hypotension, neurogenic bladder); (5) cardiac disease characterized by thickened ventricular walls in the absence of hypertension; (6) advanced atrio-ventricular block of unknown origin, particularly when accompanied by a thickened heart; and (6) vitreous body inclusions of the cotton-wool type. See Ando et al., Orphanet Journal of Rare Diseases 8:31 (2013). Definitive diagnosis of TTR amyloidosis, however, typically relies on target organ biopsies, followed by histological staining of the excised tissue with the amyloid-specific dye, Congo red. If a positive test for amyloid is observed, immunohistochemical staining for TTR is subsequently performed to ensure that the precursor protein responsible for amyloid formation is indeed TTR. For familial forms of the diseases, demonstration of a mutation in the gene encoding TTR is then needed before a definitive diagnosis can be made.

[0159] The identification of the subject can occur in a clinical setting, or elsewhere, such as in the subject's home, for example, through the subject's own use of a self-testing kit. For example, the subject can be identified based on various symptoms such as peripheral neuropathy (sensory and motor), autonomic neuropathy, gastrointestinal impairment, cardiomyopathy, nephropathy, or ocular deposition. See Ando et al., Orphanet Journal of Rare Diseases 8:31 (2013). The subject can also be identified by increased levels of non-native forms of TTR in plasma samples from the subject compared to control samples, as disclosed in the examples.

[0160] As warranted by family history, genetic testing, or medical screening for TTR amyloidosis, treatment can begin at any age (e.g., 20, 30, 40, 50, 60, or 70 years of age). Treatment typically entails multiple dosages over a period of time and can be monitored by assaying antibody or activated T-cell or B-cell responses to a therapeutic agent (e.g., a truncated form of TTR comprising amino acid residues 89-97) over time. If the response falls, a booster dosage is indicated.

[0161] In prophylactic applications, an antibody or a pharmaceutical composition of the same is administered to a subject susceptible to, or otherwise at risk of a disease (e.g., TTR amyloidosis) in a regime (dose, frequency and route of administration) effective to reduce the risk, lessen the severity, or delay the onset of at least one sign or symptom of the disease. In therapeutic applications, an antibody or immunogen to induce an antibody is administered to a subject suspected of, or already suffering from a disease (e.g., TTR amyloidosis) in a regime (dose, frequency and route of administration) effective to ameliorate or at least inhibit further deterioration of at least one sign or symptom of the disease.

[0162] A regime is considered therapeutically or prophylactically effective if an individual treated subject achieves an

outcome more favorable than the mean outcome in a control population of comparable subjects not treated by methods disclosed herein, or if a more favorable outcome is demonstrated for a regime in treated subjects versus control subjects in a controlled clinical trial (e.g., a phase II, phase II/III, or phase III trial) or an animal model at the  $p < 0.05$  or  $0.01$  or even  $0.001$  level.

5 **[0163]** An effective regime of an antibody can be used for, e.g., inhibiting or reducing aggregation of TTR in a subject having or at risk of a condition associated with TTR accumulation; inhibiting or reducing TTR fibril formation in a subject having or at risk of a condition associated with TTR accumulation; reducing or clearing TTR deposits or aggregated TTR in a subject having or at risk of a condition associated with TTR accumulation; stabilizing non-toxic conformations of TTR in a subject having or at risk of a condition associated with TTR accumulation; inhibiting toxic effects of TTR aggregates, fibrils or deposits in a subject having or at risk of a condition associated with TTR accumulation; diagnosing the presence or absence of TTR amyloid accumulation in a tissue suspected of comprising the amyloid accumulation; determining a level of TTR deposits in a subject by detecting the presence of bound antibody in the subject following administration of the antibody; detecting the presence of monomeric, misfolded, aggregated, or fibril forms of TTR in a subject; monitoring and evaluating the efficacy of therapeutic agents being used to treat patients diagnosed with a TTR amyloidosis; inducing an immune response comprising antibodies to TTR in a subject; delaying the onset of a condition associated with TTR amyloid accumulation in a subject; or treating or effecting prophylaxis of a TTR amyloidosis in a patient.

10 **[0164]** Effective doses vary depending on many different factors, such as means of administration, target site, physiological state of the subject, whether the subject is human or an animal, other medications administered, and whether treatment is prophylactic or therapeutic.

15 **[0165]** An exemplary dose range for antibodies can be from about  $0.1$ - $20$ , or  $0.5$ - $5$  mg/kg body weight (e.g.,  $0.5$ ,  $1$ ,  $2$ ,  $3$ ,  $4$  or  $5$  mg/kg) or  $10$ - $1500$  mg as a fixed dosage. The dosage depends on the condition of the patient and response to prior treatment, if any, whether the treatment is prophylactic or therapeutic and whether the disorder is acute or chronic, among other factors.

20 **[0166]** Antibody can be administered in such doses daily, on alternative days, weekly, fortnightly, monthly, quarterly, or according to any other schedule determined by empirical analysis. An exemplary treatment entails administration in multiple doses over a prolonged period, for example, of at least six months. Additional exemplary treatment regimes entail administration once per every two weeks or once a month or once every 3 to 6 months.

25 **[0167]** Antibodies can be administered via a peripheral route. Routes of administration include topical, intravenous, oral, subcutaneous, intraarterial, intracranial, intrathecal, intraperitoneal, intranasal or intramuscular. Routes for administration of antibodies can be intravenous or subcutaneous. Intravenous administration can be, for example, by infusion over a period such as  $30$ - $90$  min. This type of injection is most typically performed in the arm or leg muscles. In some methods, agents are injected directly into a particular tissue where deposits have accumulated, for example intracranial injection.

30 **[0168]** Pharmaceutical compositions for parenteral administration can be sterile and substantially isotonic ( $250$ - $350$  mOsm/kg water) and manufactured under GMP conditions. Pharmaceutical compositions can be provided in unit dose form (i.e., the dose for a single administration). Pharmaceutical compositions can be formulated using one or more physiologically acceptable carriers, diluents, excipients or auxiliaries. The formulation depends on the route of administration chosen. For injection, antibodies can be formulated in aqueous solutions, e.g., in physiologically compatible buffers such as Hank's solution, Ringer's solution, or physiological saline or acetate buffer (to reduce discomfort at the site of injection). The solution can contain formulatory agents such as suspending, stabilizing and/or dispersing agents. Alternatively antibodies can be in lyophilized form for constitution with a suitable vehicle, e.g., sterile pyrogen-free water, before use.

35 **[0169]** The regimes can be administered in combination with another agent effective in treatment or prophylaxis of the disease being treated. Such agents can include siRNA to inhibit expression of TTR or Vyndaqel, a stabilizer of TTR in tetramer formation.

40 **[0170]** After treatment, the subject's condition can be evaluated to determine the progress or efficacy of such treatment. Such methods preferably test for changes in TTR amyloid levels or levels of non-native forms of TTR. For example, TTR amyloid levels may be evaluated to determine improvement relative to the subject's TTR amyloid levels under comparable circumstances prior to treatment. The subject's TTR amyloid levels can also be compared with control populations under comparable circumstances. The control populations can be similarly afflicted, untreated subjects or normal untreated subjects (among other control subjects). Improvement relative to similarly afflicted, untreated subjects or levels approaching or reaching the levels in untreated normal subjects indicates a positive response to treatment.

45 **[0171]** TTR amyloid levels can be measured by a number of methods, including imaging techniques. Examples of suitable imaging techniques include PET scanning with radiolabeled TTR or fragments thereof, TTR antibodies or fragments thereof, Congo-red-based amyloid imaging agents, such as, e.g., PIB (US 2011/0008255), amyloid-imaging peptide p31 (Biodistribution of amyloid-imaging peptide, p31, correlates with amyloid quantitation based on Congo red tissue staining, Wall *et al.*, Abstract No. 1573, 2011 ISNM Annual Meeting), and other PET labels. Levels of non-native

forms of TTR can be measured, for example, by performing SDS-PAGE/Western blot or Meso Scale Discovery plate assays with the antibodies disclosed herein on plasma samples or biopsy samples from a subject and comparing to control samples, as described in the examples.

#### 5 **A. Diagnostics and Monitoring Methods**

10 **[0172]** Also provided are methods of detecting an immune response against TTR in a patient suffering from or susceptible to diseases associated with TTR deposition or pathogenic forms of TTR (e.g., monomeric, misfolded, aggregated, or fibril forms of TTR). The methods can be used to monitor a course of therapeutic and prophylactic treatment with the agents provided herein. The antibody profile following passive immunization typically shows an immediate peak in antibody concentration followed by an exponential decay. Without a further dose, the decay approaches pretreatment levels within a period of days to months depending on the half-life of the antibody administered. For example, the half-life of some human antibodies is of the order of 20 days.

15 **[0173]** In some methods, a baseline measurement of antibody to TTR in the subject is made before administration, a second measurement is made soon thereafter to determine the peak antibody level, and one or more further measurements are made at intervals to monitor decay of antibody levels. When the level of antibody has declined to baseline or a predetermined percentage of the peak less baseline (e.g., 50%, 25% or 10%), administration of a further dose of antibody is administered. In some methods, peak or subsequent measured levels less background are compared with reference levels previously determined to constitute a beneficial prophylactic or therapeutic treatment regime in other subjects. If the measured antibody level is significantly less than a reference level (e.g., less than the mean minus one or, preferably, two standard deviations of the reference value in a population of subjects benefiting from treatment) administration of an additional dose of antibody is indicated.

20 **[0174]** Also provided are methods of detecting monomeric, misfolded, aggregated, or fibril forms of TTR in a subject, for example, by measuring TTR amyloid or pathogenic forms of TTR (e.g., monomeric, misfolded, aggregated, or fibril forms of TTR) in a sample from a subject or by *in vivo* imaging of TTR in a subject. Such methods are useful to diagnose or confirm diagnosis of diseases associated with such pathogenic forms of TTR (e.g., TTR amyloidosis), or susceptibility thereto. The methods can also be used on asymptomatic subjects. The presence of monomeric, misfolded, aggregated, or fibril forms of TTR indicates susceptibility to future symptomatic disease. The methods are also useful for monitoring disease progression and/or response to treatment in subjects who have been previously diagnosed with a TTR amyloidosis.

25 **[0175]** Biological samples obtained from a subject having, suspected of having, or at risk of having a TTR amyloidosis can be contacted with the antibodies disclosed herein to assess the presence of monomeric, misfolded, aggregated, or fibril forms of TTR. For example, levels of monomeric, misfolded, aggregated, or fibril forms of TTR in such subjects may be compared to those present in healthy subjects. Alternatively, levels of TTR amyloid or pathogenic forms of TTR (e.g., monomeric, misfolded, aggregated, or fibril forms of TTR) in such subjects receiving treatment for the disease may be compared to those of subjects who have not been treated for a TTR amyloidosis. Some such tests involve a biopsy of tissue obtained from such subjects. ELISA assays may also be useful methods, for example, for assessing levels of monomeric, misfolded, aggregated, or fibril forms of TTR in fluid samples. Some such ELISA assays involve anti-TTR antibodies that preferentially bind monomeric, misfolded, aggregated, or fibril forms of TTR relative to normal tetrameric forms of TTR.

30 **[0176]** The *in vivo* imaging methods can work by administering a reagent, such as antibody that binds to monomeric, misfolded, aggregated, or fibril forms of TTR in the subject, and then detecting the reagent after it has bound. Such antibodies typically bind to an epitope within residues 89-97 of TTR. If desired, the clearing response can be avoided by using antibody fragments lacking a full length constant region, such as Fabs. In some methods, the same antibody can serve as both a treatment and diagnostic reagent.

35 **[0177]** Diagnostic reagents can be administered by intravenous injection into the body of the subject, or via other routes deemed reasonable. The dose of reagent should be within the same ranges as for treatment methods. Typically, the reagent is labeled, although in some methods, the primary reagent with affinity for monomeric, misfolded, aggregated, or fibril forms of TTR is unlabeled and a secondary labeling agent is used to bind to the primary reagent. The choice of label depends on the means of detection. For example, a fluorescent label is suitable for optical detection. Use of paramagnetic labels is suitable for tomographic detection without surgical intervention. Radioactive labels can also be detected using PET or SPECT.

40 **[0178]** Diagnosis is performed by comparing the number, size, and/or intensity of labeled loci to corresponding base line values. The base line values can represent the mean levels in a population of undiseased individuals. Base line values can also represent previous levels determined in the same subject. For example, base line values can be determined in a subject before beginning treatment, and measured values thereafter compared with the base line values. A decrease in values relative to base line generally signals a positive response to treatment.

## IX. Kits

**[0179]** The disclosure further provides kits (e.g., containers) comprising an antibody disclosed herein and related materials, such as instructions for use (e.g., package insert). The instructions for use may contain, for example, instructions for administration of the antibody and optionally one or more additional agents. The containers of antibody may be unit doses, bulk packages (e.g., multi-dose packages), or sub-unit doses.

**[0180]** Package insert refers to instructions customarily included in commercial packages of therapeutic products that contain information about the indications, usage, dosage, administration, contraindications and/or warnings concerning the use of such therapeutic products

**[0181]** Kits can also include a second container comprising a pharmaceutically-acceptable buffer, such as bacteriostatic water for injection (BWFJ), phosphate-buffered saline, Ringer's solution and dextrose solution. It can also include other materials desirable from a commercial and user standpoint, including other buffers, diluents, filters, needles, and syringes.

## X. Other Applications

**[0182]** The antibodies can be used for detecting monomeric, misfolded, aggregated, or fibril forms of transthyretin (TTR), or fragments thereof, in the context of clinical diagnosis or treatment or in research. For example, the antibodies can be used to detect the presence of monomeric, misfolded, aggregated, or fibril forms of TTR in a biological sample as an indication that the biological sample comprises TTR amyloid deposits. Binding of the antibodies to the biological sample can be compared to binding of the antibodies to a control sample. The control sample and the biological sample can comprise cells of the same tissue origin. Control samples and biological samples can be obtained from the same individual or different individuals and on the same occasion or on different occasions. If desired, multiple biological samples and multiple control samples are evaluated on multiple occasions to protect against random variation independent of the differences between the samples. A direct comparison can then be made between the biological sample(s) and the control sample(s) to determine whether antibody binding (*i.e.*, the presence of monomeric, misfolded, aggregated, or fibril forms of TTR) to the biological sample(s) is increased, decreased, or the same relative to antibody binding to the control sample(s). Increased binding of the antibody to the biological sample(s) relative to the control sample(s) indicates the presence of monomeric, misfolded, aggregated, or fibril forms of TTR in the biological sample(s). In some instances, the increased antibody binding is statistically significant. Optionally, antibody binding to the biological sample is at least 1.5-fold, 2-fold, 3-fold, 4-fold, 5-fold, 10-fold, 20-fold, or 100-fold higher than antibody binding to the control sample.

**[0183]** In addition, the antibodies can be used to detect the presence of monomeric, misfolded, aggregated, or fibril forms of TTR in a biological sample to monitor and evaluate the efficacy of a therapeutic agent being used to treat a patient diagnosed with a TTR amyloidosis. A biological sample from a patient diagnosed with a TTR amyloidosis is evaluated to establish a baseline for the binding of the antibodies to the sample (*i.e.*, a baseline for the presence of the monomeric, misfolded, aggregated, or fibril forms of TTR in the sample) before commencing therapy with the therapeutic agent. In some instances, multiple biological samples from the patient are evaluated on multiple occasions to establish both a baseline and measure of random variation independent of treatment. A therapeutic agent is then administered in a regime. The regime may include multiple administrations of the agent over a period of time. Optionally, binding of the antibodies (*i.e.*, presence of monomeric, misfolded, aggregated, or fibril forms of TTR) is evaluated on multiple occasions in multiple biological samples from the patient, both to establish a measure of random variation and to show a trend in response to immunotherapy. The various assessments of antibody binding to the biological samples are then compared. If only two assessments are made, a direct comparison can be made between the two assessments to determine whether antibody binding (*i.e.*, presence of monomeric, misfolded, aggregated, or fibril forms of TTR) has increased, decreased, or remained the same between the two assessments. If more than two measurements are made, the measurements can be analyzed as a time course starting before treatment with the therapeutic agent and proceeding through the course of therapy. In patients for whom antibody binding to biological samples has decreased (*i.e.*, the presence of monomeric, misfolded, aggregated, or fibril forms of TTR), it can be concluded that the therapeutic agent was effective in treating the TTR amyloidosis in the patient. The decrease in antibody binding can be statistically significant. Optionally, binding decreases by at least 1%, 2%, 3%, 4%, 5%, 10%, 15%, 20%, 30%, 40%, 50%, 60%, 70%, 80%, 90%, or 100%. Assessment of antibody binding can be made in conjunction with assessing other signs and symptoms of TTR amyloidosis.

**[0184]** The antibodies can also be used as research reagents for laboratory research in detecting monomeric, misfolded, aggregated, or fibril forms of TTR, or fragments thereof. In such uses, antibodies can be labeled with fluorescent molecules, spin-labeled molecules, enzymes, or radioisotopes, and can be provided in the form of kit with all the necessary reagents to perform the detection assay. The antibodies can also be used to purify monomeric, misfolded, aggregated, or fibril forms of TTR, or binding partners of monomeric, misfolded, aggregated, or fibril forms of TTR, e.g., by affinity chromatography.

**[0185]** The antibodies can also be used for inhibiting or reducing aggregation of TTR, inhibiting or reducing TTR fibril formation, reducing or clearing TTR deposits or TTR aggregates, or stabilizing non-toxic conformations of TTR in a biological sample. The biological sample can comprise, for example, blood, serum, plasma, or tissue (e.g., tissue from the heart, peripheral nervous system, autonomic nervous system, kidneys, eyes, or gastrointestinal tract). In some instances, TTR aggregation, TTR fibril formation, or TTR deposits are inhibited or reduced by at least 10%, 20%, 25%, 30%, 40%, 50%, or 75%, (e.g., 10%-75% or 30%-70%). Assays for detecting fibril formation are described elsewhere herein. See also US 2014/0056904.

**[0186]** If different versions of a sequence are associated with an accession number at different times, the version associated with the accession number at the effective filing date of this application is meant. The effective filing date means the earlier of the actual filing date or filing date of a priority application referring to the accession number if applicable. Likewise if different versions of a publication, website or the like are published at different times, the version most recently published at the effective filing date of the application is meant unless otherwise indicated. Any feature, step, element, embodiment, or aspect of the disclosure can be used in combination with any other unless specifically indicated otherwise. Although the present invention has been described in some detail by way of illustration and example for purposes of clarity and understanding, it will be apparent that certain changes and modifications may be practiced within the scope of the appended claims.

## EXAMPLES

### Example 1. Identification of mis-TTR Monoclonal Antibodies

**[0187]** Conformationally-specific monoclonal antibodies against monomeric, mis-folded, fibril, or aggregated forms of TTR (mis-TTR) were generated, screened, expressed, and purified as described in Materials and Methods (a-d). In order to generate mis-TTR monoclonal antibodies, the crystal structure of human tetrameric TTR was examined to find regions of the protein that are buried in the tetramer, but exposed upon dissociation of the tetramer into its monomeric subunits. The region identified was residues 89-97 (EHA<sub>2</sub>EVVFTA) (SEQ ID NO: 113) located within the F strand of TTR and sequestered at the dimer interface of the tetrameric protein. A BLAST search of the protein database did not reveal any other human proteins possessing this sequence.

**[0188]** A peptide comprising this sequence (ggEHA<sub>2</sub>EWFTA<sub>2</sub>ggkg) (SEQ ID NO: 114), was synthesized. Capitalized letters represent residues 89-97 of TTR. Lower case letters represent additional linker residues added to increase the solubility of the antigenic peptide and to establish the 9 amino acid fragment as an internal sequence. This peptide was linked to a poly-lysine dendritic core, generating a multiple antigenic peptide immunogen (TTR-MAP) comprising a core of lysine residues with multiple branches linked to the TTR 89-97 peptide. The antibodies listed in Table 2 were generated against TTR-MAP

**[0189]** In addition to this multiple antigenic peptide, two other immunogens containing the same TTR fragment were generated by covalently linking similar TTR 89-97 peptides (AccggEHA<sub>2</sub>EVVFTA-amide (SEQ ID NO: 115) and Ac-EHA<sub>2</sub>EWFTA<sub>2</sub>ggg-amide (SEQ ID NO: 116)) via the N- and C-terminal cysteine residues to keyhole limpet hemocyanin (TTR89-97-N-KLH and TTR89-97-C-KLH).

**[0190]** Following antibody generation, screening, expression, and purification, detailed binding kinetic parameters (association rate ( $k_a$ ), dissociation rate ( $k_d$ ), and binding affinity constant ( $K_D$ )) were determined for lead mis-TTR antibodies by Surface Plasmon Resonance (SPR) for recombinant human TTR F87M/L110M, as shown in Table 2. Anti-mouse IgG (GE Healthcare) was immobilized on a sensor chip C5 (lacking dextran chains) via amine coupling following the instructions provided in the GE Healthcare anti-mouse kit, and mis-TTR mAbs were captured to a level to ensure a maximum binding of analyte of 30-50 RU. Various concentrations of analyte (recombinant human TTR F87M/L110M) were passed over the captured ligand at 30  $\mu$ l/min in running buffer (HBS + 0.05% P-20, 1 mg/mL BSA) in 3-fold dilutions. For each concentration, the reaction proceeded for a time frame allowing for the higher analyte concentrations to reach equilibrium during association, as well as at least 10% of signal to decay during dissociation. At least one concentration (not the highest or lowest) was run in duplicate. Concentration ranges of analyte were selected based on preliminary experimentation to span at least 10-fold above  $K_D$  to 10-fold below  $K_D$ .

**[0191]** The results of SPR analysis of lead mis-TTR mAbs is shown in Table 2 below.

**Table 2**

<b>SPR Analysis of Lead mis-TTR Antibodies Binding to Human TTR (F87M/L110M)</b>				
mAb	$k_a$ (1/Ms)	$k_d$ (1/s)	$K_D$ (M)	Rmax
9D5	2.715E+4	4.930E-4	1.816E-8	31.55
14G8	2.880E+4	5.358E-4	1.861E-8	27.13

(continued)

SPR Analysis of Lead mis-TTR Antibodies Binding to Human TTR (F87M/L110M)				
mAb	$k_a$ (1/Ms)	$k_d$ (1/s)	$K_D$ (M)	Rmax
5A1	6.107E+4	4.693E-4	7.684E-9	30.98
6C1	4.607E+4	4.151E-4	9.010E-9	26.32

### Example 2. Binding of mis-TTR Antibodies to TTR Antigen

**[0192]** Four lead mis-TTR mAbs (9D5, 14G8, 6C1, and 5A1) were assayed by ELISA at concentrations ranging from 0.31 to 2.5  $\mu\text{g/ml}$  using both pH4.0-treated TTR (pH4-TTR) and native TTR as the coating antigen. TTR antigen preparation and ELISA protocols are described elsewhere in Materials and Methods (e-g).

**[0193]** The resulting binding curves and tabulated  $K_d$  and  $B_{\text{max}}$  values are shown in Figure 3 and Table 3 below. The results in Figure 3 are presented in arbitrary units (a.u.) on the y-axis. All mAbs showed significant binding to pH4-TTR with  $K_d$  values ranging from 16 nM (6C1) to 282 nM (9D5).  $B_{\text{max}}$  values for binding to pH4-TTR ranged from a low of 0.65 a.u. (14G8) to a high of 2.02 (9D5). In contrast to the binding to pH4-TTR, none of the antibodies showed significant binding to native TTR, indicating that all TTR antibodies generated were specific for non-native forms of TTR.

Table 3

ELISA Analysis of Lead mis-TTR Antibodies Binding to pH4-TTR		
mAb	$K_d$ (nM)	$B_{\text{max}}$ (a.u.)
9D5	282	2.02
14G8	108	0.65
6C1	16	1.07
5A1	23	1.61

### Example 3. Analysis of mis-TTR Antibodies by SDS-PAGE and Native-PAGE

**[0194]** 9D5 and 14G8 were analyzed by SDS-PAGE/Western to demonstrate specificity of binding toward monomeric/denatured forms of TTR versus native, non-denatured TTR. SDS-PAGE, Native-PAGE, and Western Blot protocols are described elsewhere in the Methods and Materials (h-j).

**[0195]** Non-denatured TTR or pH4-TTR was run on an SDS-PAGE gel alongside heat-denatured TTR and heat-denatured pH4-TTR. After electrophoresis, the gel was Western blotted onto nitrocellulose and stained with TTR mAbs 9D5 and 14G8. Both antibodies only recognized TTR when it was treated at pH4 or when TTR or pH4-TTR was first heat-denatured prior to SDS-PAGE. These 9D5 and 14G8 thus show a specificity for TTR conformers generated either by denaturation of TTR or by treatment of TTR at pH4.

**[0196]** 6C1 and 5A1 along with total TTR mAbs (7G7, 8C3) and the commercially available Sigma polyclonal antibody were also analyzed by SDS-PAGE/Western. Each blot contained stained molecular weight markers, non-denatured TTR, and pH4-TTR.

**[0197]** The stained SDS-PAGE gel showed that the major species present in the non-denatured TTR sample was an  $\sim 38$  kDa dimer. In contrast, the major component present in the pH4-TTR sample ran as an  $\sim 35$  kDa dimer with a small amount of dimer of an  $\sim 15$  kDa monomer. This dimer ran as a slightly smaller protein than the dimer present in the non-denatured TTR sample, indicating a conformational difference between these two TTR dimer species.

**[0198]** The Western blots of TTR and pH4-TTR using the four mis-TTR antibodies showed that these mAbs do not recognize non-denatured TTR, but do bind to both the denatured monomer and dimer present in the pH4-TTR sample. Thus, the four mis-TTR mAbs (9D5, 14G8, 6C1, and 5A1) show similar specificities for non-native conformations of TTR when analyzed by SDS-PAGE/Westerns.

**[0199]** In contrast to the four mis-TTR mAbs, the two TTR control mAbs, 7G7 and 8C3, generated through immunization of mice with intact TTR recognized all TTR species present in the TTR and pH4-TTR samples, including tetrameric TTR species. Thus unlike the mis-TTR mAbs, these control mAbs bind TTR but with no conformational specificity. The Sigma polyclonal antibody behaved similarly to the 7G7 and 8C3 control mAbs.

**[0200]** TTR and pH4-TTR were also run on a native gel to see if the four mis-TTR mAbs were capable of showing conformation specificity under non-denaturing gel conditions. On a stained native PAGE gel, TTR ran as an  $\sim 35$  kDa

native dimer with a small amount of tetramer. In contrast, pH4-TTR ran primarily as a high molecular-weight smear with a trace amount of the ~35 kDa dimer. The non-specific Sigma polyclonal antibody recognized all TTR species present in both the TTR and the pH4-TTR sample. In contrast, 9D5 only recognized the high molecular weight TTR species present in the pH4-TTR sample. As observed in the SDS-PAGE/Western study, 9D5 did not recognize any of the native TTR species.

[0201] All four mis-TTR mAbs were subsequently analyzed by native-PAGE/Western blot. As expected and similar to 9D5, the other mis-TTR mAbs, 14G8, 6C1, and 5A1, specifically bound to the high molecular weight non-native forms of TTR present in the pH4-TTR sample. None of these antibodies recognized the ~35 kDa native TTR dimer. These results indicate that the four mis-TTR mAbs behave similarly and recognize only non-native TTR species that are conformationally distinct from native TTR.

#### Example 4. Inhibition of TTR Fiber Formation by mis-TTR Antibodies

[0202] TTR-Y78F is a TTR variant containing a point mutation at position 78 in the protein sequence that destabilizes the TTR tetramer. With time and under mildly acidic conditions, this TTR variant dissociates into its monomeric subunits which can then go on to aggregate and form fibers capable of binding to thioflavin-T. The extent of fiber formation can thus be monitored by measuring thioflavin-T fluorescence at 480nm. Introduction of a mis-TTR antibody specific for dissociated TTR monomers or aggregates would prevent the assembly of TTR fibers resulting in a decrease in thioflavin-T fluorescence relative to a no-antibody control reaction. Protocols for examining inhibition of TTR fiber formation are described elsewhere in the Materials and Methods (k).

[0203] All four mis-TTR antibodies strongly inhibited the formation of thioflavin-T reactive TTR-Y78F fibers relative to the isotype control. The results are shown in Figure 4A and are presented in arbitrary units (a.u.) on the y-axis. Mis-TTR antibody 5A1 almost completely inhibited fiber formation. These results are consistent with the notion that mis-TTR antibodies bind monomeric and/or aggregated forms of TTR, thereby preventing the formation of TTR fibers.

[0204] Table 4 summarizes the characterization data obtained for the set of 4 mis-TTR antibodies (9D5, 14G8, 6C1, and 5A1) that showed good conformational selectivity for non-native forms of TTR. These antibodies had affinities ( $K_D$ ) for pH4-TTR ranging from 14.5 nM (6C1) to 257 nM (9D5) and  $B_{max}$  values ranging from 0.65 a.u. (14G8) to 2.02 (9D5). None of these antibodies recognized native TTR, but did bind to pH4-TTR on an SDS-PAGE/Western and to the high molecular weight TTR aggregates on a native-PAGE/Western. These antibodies also inhibited the formation of TTR fibrils in the fibril formation assay using Thio-T as the readout.

Table 4

mis-TTR-Y78F mAb Characterization Summary Table						
Clone ID	Sandwich ELISA (pH4-TTR)		Western Blot			% Inh. Fibrils (Thio- T)
	$K_D$ (nM)	B max (OD <sub>450</sub> a.u.)	SDS-PAGE		Native (HMW-TTR)	
			(TTR)	(pH4- TTR)		
9D5	257	2.02	-	+++	+++	83
14G8	98.7	0.65	-	+++	++	65
6C1	14.6	1.07	-	+++	+++	72
5A1	21.3	1.61	-	+++	+++	100

[0205] TTR-V122I is a TTR variant containing a single point mutation at position 122 that destabilizes the tetramer. Fibril formation was associated with an increase in ThT fluorescence. Increasing 14G8 mAb concentrations caused a monotonic decrease in ThT fluorescence indicating a substoichiometric inhibition of TTR fibrillation ( $IC_{50}=0.028\pm 0.009$  mg/mL; n=3; Figure 4B and Table 4a). The isotype control mAb did not cause inhibition of TTR fibrillation (Figure 4C), thus demonstrating the specificity of 14G8 mediated inhibition.

[0206] Comparable substoichiometric  $IC_{50}$  values determined for 5A1 and 6C1 (Table 4a) suggested analogous mechanisms of fibril inhibition for each of these mis-TTR mAbs. In contrast, 9D5 unexpectedly failed to inhibit TTR-V122I fibril formation, despite showing similar specificity and affinity for non-native TTR. It remains to be explored whether 9D5 is more sensitive to the assay conditions used.

**Table 4a**  
**mis-TTR-V122I mAb Characterization Summary**  
**Table**

Antibody	IC <sub>50</sub> ± SD (mg/mL)
9D5	No inhibition
14G8	0.028 ± 0.009
6C1	0.048 ± 0.059
5A1	0.015 ± 0.02
EG 27/1	No inhibition

**Example 5. Immunohistochemical (IHC) Characterization of ATTR Tissue Using mis-TTR mAbs**

**[0207]** The lead mis-TTR mAbs raised to the TTR 89-97 fragment of the transthyretin protein were immunohistochemically tested on fresh frozen and paraffin processed tissue from confirmed TTR cardiac amyloidosis patients. Protocols for obtaining and preparing cardiac tissue samples, immunohistochemistry (IHC), and image analysis, are provided elsewhere in the Materials and Methods (I-o). The antibodies used for IHC are described in Table 5.

**Table 5**

<b>Antibodies Used for Immunohistochemical Characterization</b>				
<b>Antibody</b>	<b>Antibody Type</b>	<b>Vendor</b>	<b>Stain Cardiac Tissue</b>	<b>Concentration</b>
14G8	mis-TTR	Prothena Biosciences	Yes	0.5 µg/mL
9D5	mis-TTR	Prothena Biosciences	Yes	0.5 µg/mL
6C1	mis-TTR	Prothena Biosciences	Yes	0.5 µg/mL
5A1	mis-TTR	Prothena Biosciences	Yes	0.5 µg/mL
7G7	TTR	Prothena Biosciences	Yes	0.5 µg/mL
6F10	Isotype Control	Prothena Biosciences	No	0.5 µg/mL
Prealbumin (A0002)	TTR	Dako North America	Yes	1:2,000 & 1:20,000
Kappa Light Chains (A0191)	LC-κ	Dako North America	No	1:8,000
Lambda Light Chains (A0193)	LC-λ	Dako North America	No	1:8,000
Amyloid A (M0759)	AA	Dako North America	No	1:8,000

**[0208]** Cardiac tissue samples were obtained from patients with confirmed diagnoses of ATTR mutations. Demographics for cases examined immunohistochemically were as follows and are summarized in Table 6: FAC = familial amyloidotic cardiomyopathy; FAP = familial amyloidotic polyneuropathy; 1° AL = light-chain amyloidosis; ATTR = transthyretin-mediated amyloidosis; Unk = Unknown

Table 6

Immunohistochemical Staining of Cardiac Tissue Samples with mis-TTR Antibodies				
Case	Diagnosis	TTR Mutations	Format	Stained with TTR Antibodies?
Patient 1	FAC	Ileu122	Frozen	Yes
Patient 2	FAP	Wild type	Frozen	Yes
Patient 3	FAP	84Ser	Frozen	Yes
Patient 4	FAP	84Ser	Frozen	Yes
Patient 5	1° AL	-	Frozen	No
Patient 6	1° AL	-	Frozen	No
Patient 7	ATTR	10Arg	Frozen	Yes
Patient 8	ATTR	V122I	Frozen	Yes
Patient H1	ATTR	Val122Ile	FFPE	Yes
Patient H2	ATTR	Thr60Ala	FFPE	Yes
Patient H3	ATTR	Thr49Ala	FFPE	Yes
Patient H4	ATTR	Ile84Ser	FFPE	Yes
Patient H5	Unk.	Senile Cardiac	FFPE	Yes
Patient H6	ATTR	Ile84Ser	FFPE	Yes

**[0209]** Mouse monoclonal antibodies (mis-TTR mAbs) raised to the 89-97 fragment of the transthyretin protein were immunohistochemically tested on fresh frozen and paraffin processed tissue from confirmed TTR cardiac amyloidosis patients. Each mis-TTR antibody showed strong immunoreactivity on ATTR cardiac tissue. Dark staining was observed in deposits throughout the myocardium and the vasculature. When immunoreactivity was compared to staining with Congo Red or Thioflavin-T, the majority of the immunoreactivity in the tissue showed high congruence with Congo red birefringence and Thioflavin T-positive staining. This confirms the beta pleated sheet nature of the TTR amyloid deposited in this tissue. 14G8, 9D5, 6C1, and 5A1 also detected pre-amyloid TTR present in areas of the myocardium that were TTR-immunopositive but Congo red or Thioflavin T-negative. Both the IgG-isotype control antibody and primary antibody omission sections were negative for staining across all tissues tested. Antibodies reactive toward other amyloidogenic proteins (lambda and kappa light chains or amyloid A) were non-reactive toward the ATTR cardiac tissue used in this analysis, indicating that deposits were specifically TTR in nature.

**[0210]** The staining patterns of mis-TTR antibodies were compared to that obtained with a well-characterized commercial TTR reference antibody (prealbumin, A0002; Dako; Carpinteria, CA). The DAKO reference antibody stained the diseased myocardium in the same areas as the mis-TTR antibodies, but produced a more diffuse staining pattern. The DAKO reference antibody did not stain the congophillic TTR amyloid deposits present on the vasculature as strongly as the mis-TTR antibodies.

**[0211]** The mis-TTR antibodies did not stain normal, non-disease tissue. Furthermore, as expected, staining with an isotype control antibody, 6F10, was also negative.

**[0212]** In order to determine if the reactivity of mis-TTR antibodies was specific for TTR deposits, cross reactivity of these antibodies toward cardiac tissue derived from patients diagnosed with primary AL amyloidosis was examined. As expected, no staining of AL amyloid tissue was observed, confirming that TTR antibodies react specifically toward ATTR diseased tissue.

**[0213]** Cardiac tissue from patients with confirmed diagnoses of senile systemic amyloidosis or from patients with confirmed FAC, or FAP caused by point mutations in the TTR gene also stained positively with 14G8, 9D5, 6C1, and 5A1. These results indicate that mis-TTR antibodies have the ability to recognize TTR deposits in cardiac tissue regardless of the ATTR genotype.

**[0214]** Other non-cardiac tissues known to express TTR were also examined for staining by mis-TTR antibodies and compared to the staining obtained using the DAKO reference antibody. As expected, the liver, pancreas and choroid plexus all stained positively for TTR using the Dako reference antibody. In contrast, 14G8 only stained the pancreatic alpha cells located in the islets of Langerhans and the choroid plexus, suggesting that some of the TTR localized to these organs are conformationally distinct from TTR expressed in the liver. The lack of mis-TTR mAb immunoreactivity in the liver suggests that the large amount of TTR expressed there is primarily tetrameric, native TTR and does not have

the exposed mis-TTR epitope. Similar results were obtained when the same tissues were stained with 9D5, 6C1, and 5A1.

#### Example 6. Analysis of ATTR vs Normal Human Plasma by SDS-PAGE/Western Blot and by Meso Scale Discovery (MSD) Plate Assay

[0215] Six plasma samples from patients confirmed for V30M ATTR (#21, #22, #23, #24, #25, #27) and 6 samples from normal subjects (Sample #11, #12, #15, #18, #19, #20) were obtained from M. Saraiva (Porto University, Portugal). Sample #C6 was a normal human serum sample obtained from a commercial source (Bioreclamation/VT). Samples were analyzed by SDS-PAGE and Western blot, or by MesoScale Discovery (MSD) Plate Assay. Protocols for these assays are described elsewhere in the Materials and Methods (p-r). A standard curve was generated for the MSD Plate Assay using 6C1.

[0216] In the resulting Western blots using the 9D5 and the 5A1 mis-TTR mAbs, differences in banding patterns between normal and TTR-V30M diseased plasma samples could be detected. All plasma samples contained an ~14 kDa TTR band that co-migrated with the non-native TTR monomer present in the pH4-TTR reference sample. In general, plasma samples derived from TTR-V30M patients (#21, 22, 23, 24, 25, & 27) had more of this monomeric mis-TTR species. In addition, plasma samples derived from V30M patients also contained an ~30kDa band that co-migrated with the non-native TTR dimer present in the pH4-TTR reference sample. With the exception of samples #12 and #18, plasma samples derived from normal individuals possessed less of this dimer species.

[0217] The above Western blots were scanned and the intensities of the combined 9D5- or 5A1-reactive TTR dimer and monomer bands were plotted for each sample. The results are shown in Figure 5A (9D5) and 5B (5A1) and are presented in arbitrary units (a.u.) on the y-axis. With the exception of plasma samples #15 and #18, plasma samples derived from normal individuals (11, 12, 19, and 20) contained less 9D5- and 5A1-reactive dimer and monomer species than samples derived from V30M patients (21-25 and 27).

[0218] The 12 serum samples analyzed by 9D5 and 5A1 Western blot were also analyzed by MSD plate assay using 6C1 as the mis-TTR capture antibody and the Dako-SulfoTag antibody as the detection antibody. Results of these MSD assays are shown in Figure 6 and are presented in arbitrary units (a.u.) on the y-axis. Samples 11, 12, 15, 18, 19, and 20 represent normal plasma. Samples 21-25 and 27 represent V30M diseased plasma.

[0219] With the exception of plasma samples #15 and #18, the amount of 6C1-reactive TTR present in plasma samples derived from normal individuals was lower than that observed in plasma from TTR-V30M diseased individuals. The levels of 6C1 reactivity measured by MSD assay correlated very well with the amount of 9D5 reactive dimer and monomer observed above by SDS-PAGE/Western.

[0220] In order to determine the concentration of the reactive TTR species present in plasma samples, the same samples were re-assayed using 6C1 as the capture antibody and 8C3-SulfoTag as the detection antibody. MSD signals were converted to ng/ml concentrations of reactive TTR species using the TTR F87M/L110M standard curve generated above. Based on this analysis, the average concentration of 6C1-reactive TTR present in the control samples was 271 +/- 185 ng/ml. In contrast, the average concentration of reactive TTR present in the V30M diseased plasma samples was higher, at 331 +/- 95 ng/ml. Taken together, these MSD results suggest that mis-TTR antibodies are capable of distinguishing between ATTR disease versus normal plasma. This warrants further development of mis-TTR antibodies for use in diagnostic assays of ATTR disease.

#### Example 7. Design of Humanized 9D5 Antibodies

[0221] The starting point or donor antibody for humanization was the mouse antibody 9D5. The heavy chain variable amino acid sequence of mature m9D5 is provided as SEQ ID NO: 1. The light chain variable amino acid sequence of mature m9D5 is provided as SEQ ID NO: 16. The heavy chain CDR1, CDR2, and CDR3 amino acid sequences are provided as SEQ ID NOS: 13-15, respectively (as defined by Kabat). A composite Chothia-Kabat CDR-H1 is provided as SEQ ID NO: 117. The light chain CDR1, CDR2, and CDR3 amino acid sequences are provided as SEQ ID NOS: 24-26, respectively (as defined by Kabat). Kabat numbering is used throughout in this Example.

[0222] The variable kappa (Vk) of m9D5 belongs to mouse Kabat subgroup 2, which corresponds to human Kabat subgroup 3. The variable heavy (Vh) of m9D5 belongs to mouse Kabat subgroup 3d, which corresponds to Kabat subgroup 3. See Kabat et al. Sequences of Proteins of Immunological Interest, Fifth Edition. NIH Publication No. 91-3242, 1991. The 16-residue CDR-L1 belongs to canonical class 4, the 7-residue CDR-L2 belongs to canonical class 1, and the 9-residue CDR-L3 belongs to canonical class 1 in Vk. See Martin & Thornton, J. Mol. Biol. 263:800-15, 1996. The 10-residue composite Chothia-Kabat CDR-H1 belongs to canonical class 1, and the 17-residue CDR-H2 belongs to canonical class 1. See Martin & Thornton, J Mol. Biol. 263:800-15, 1996. The CDR-H3 has no canonical classes.

[0223] The residues at the interface between the Vk and Vh domains are the ones commonly found, except that Leu is at position 47 in heavy chain, whereas Tyr typically is that this position. This position is a candidate for backmutation.

[0224] A search was made over the protein sequences in the PDB database (Deshpande et al., Nucleic Acids Res.

33: D233-7, 2005) to find structures which would provide a rough structural model of 9D5. The crystal structure of antibody fab (pdb code 1MJU) (Ruzheinikov et al., J. Mol. Biol. 332(2):423-435, 2003) was used for the Vk structure since it had good resolution (1.22 Å) and overall sequence similarity to 9D5 Vk, retaining the same canonical structure for the loop as 9D5. A monomeric antibody with pdb code 1SEQ (Covaceuszach et al., Acta Crystallogr. D Biol. Crystallogr. 57 (PT 9), 1307-1309, 2001) was used for the Vh structure since it had good sequence similarity and resolution (2.0 Å), and it has the same canonical structures for CDR-H1 and CDR-H2 as that of 9D5 VH. We modeled 9D5 Vh chain using the 1MQK structure as well, since it has a better resolution of 1.28 Å (Ostermeier et al., Proteins 21(1):74-77, 1995). BioLuminate software (licensed from Schrodinger Inc.) was used to model a rough structure of 9D5.

[0225] A search of the non-redundant protein sequence database from NCBI allowed selection of suitable human frameworks into which to graft the murine CDRs. For Vh, human Ig heavy chain BAC02114 (GI: 21670209) (SEQ ID NO: 3) was chosen (Akahori et al., Construction and characterization of antibody libraries: isolation of therapeutic human antibodies and application to functional genomics, Direct Submission, July 25, 2001). It shares the canonical forms of 9D5. It is a member of Kabat human heavy subgroup 1. 9D5 Vh has some unique framework residue and any human framework acceptor did not show very high homology. Therefore, we used a second Framework, AAX82494 (GI: 62421461) (SEQ ID NO: 4) (Lundquist et al., Infect. Immun. 74(6), 3222-3231, 2006) as well to make a hybrid acceptor framework. For Vk, a human kappa light chain with NCBI accession code ABC66952 (GI: 84798006) (Shriner et al., Vaccine 24(49-50):7159-7166, 2006) was chosen (SEQ ID NO: 18). It has the same canonical classes for CDR-L1 and L2 as that for the parental Vk. It is a member of Kabat human kappa subgroup 2.

[0226] Eight humanized heavy chain variable region variants and five humanized light chain variable region variants were constructed containing different permutations of substitutions (Hu9D5VHv1, 2, 2b, 3, 3b, 4, 4b, and 5 (SEQ ID NOS: 5-12, respectively) and Hu9D5VLv1-5 (SEQ ID NOS: 19-23, respectively)) (Tables 7 and 8). The exemplary humanized Vh and Vk designs, with backmutations and other mutations based on selected human frameworks, are shown in Tables 7 and 8, respectively. The gray-shaded areas in the first column in Tables 7 and 8 indicate the CDRs as defined by Chothia, and the gray-shaded areas in the remaining columns in Tables 7 and 8 indicate the CDRs as defined by Kabat. SEQ ID NOS: 5-12 and 19-23 contain backmutations and other mutations as shown in Table 9. The amino acids at positions H42, H47, H69, H82, H82(b), H108, L8, L9, L18, L19, L36, L39, L60, L70, and L74 in Hu9D5VHv1, 2, 2b, 3, 3b, 4, 4b, and 5 and Hu9D5VLv1-5 are listed in Tables 10 and 11.

Table 7  
Humanized 9D5 Vh Regions

Chothia residue #	Kabat residue #	Linear residue #	FR or CDR	Murine 9D5 VH (SEQ ID NO: 1)	Human VH Acceptor Fr Acc#BAC02114 (SEQ ID NO: 3)	Human VH Acceptor Fr Acc#AAX82494 (SEQ ID NO: 4)	Hu9D5 VHv1 (SEQ ID NO: 5)	Hu9D5 VHv2 (SEQ ID NO: 6)	Hu9D5 VHv3b (SEQ ID NO: 7)	Hu9D5 VHv3 (SEQ ID NO: 8)	Hu9D5 VHv3b (SEQ ID NO: 9)	Hu9D5 VHv4 (SEQ ID NO: 10)	Hu9D5 VHv4b (SEQ ID NO: 11)	Hu9D5 VHv5 (SEQ ID NO: 12)
1	1	1	FR	E	E	Q	E	E	E	E	E	E	E	E
2	2	2	FR	V	V	V	V	V	V	V	V	V	V	V
3	3	3	FR	K	Q	Q	Q	Q	Q	Q	Q	Q	Q	Q
4	4	4	FR	E	L	L	L	L	L	L	L	L	L	L
5	5	5	FR	V	V	Q	V	V	V	V	V	V	V	V
6	6	6	FR	E	E	E	E	E	E	E	E	E	E	E

**Table 7**  
**Humanized 9D5 Vh Regions**

Chothia residue #	Kabat residue #	Linear residue #	FR or CDR	Marine 9D5 VH (SEQ ID NO: 1)	Hu VH Acceptor #1 Ace#BAC03114 (SEQ ID NO: 3)	Hu VH Acceptor #r Ace#AA82494 (SEQ ID NO: 4)	Hu9D5 VHv1 (SEQ ID NO: 5)	Hu9D5 VHv2 (SEQ ID NO: 6)	Hu9D5 VHv2b (SEQ ID NO: 7)	Hu9D5 VHv3 (SEQ ID NO: 8)	Hu9D5 VHv3b (SEQ ID NO: 9)	Hu9D5 VHv4 (SEQ ID NO: 10)	Hu9D5 VHv4b (SEQ ID NO: 11)	Hu9D5 VHv5 (SEQ ID NO: 12)
7	7	7	FR	S	S	S	S	S	S	S	S	S	S	S
8	8	8	FR	G	G	G	G	G	G	G	G	G	G	G
9	9	9	FR	G	G	G	G	G	G	G	G	G	G	G
10	10	10	FR	G	G	G	G	G	G	G	G	G	G	G
11	11	11	FR	L	L	L	L	L	L	L	L	L	L	L
12	12	12	FR	V	V	V	V	V	V	V	V	V	V	V
13	13	13	FR	Q	Q	K	Q	Q	Q	Q	Q	Q	Q	Q
14	14	14	FR	F	F	P	F	F	F	F	F	F	F	F
15	15	15	FR	G	G	G	G	G	G	G	G	G	G	G
16	16	16	FR	G	G	G	G	G	G	G	G	G	G	G
17	17	17	FR	S	S	S	S	S	S	S	S	S	S	S
18	18	18	FR	L	L	L	L	L	L	L	L	L	L	L
19	19	19	FR	K	R	K	R	R	K	R	K	R	K	R
20	20	20	FR	L	L	L	L	L	L	L	L	L	L	L
21	21	21	FR	S	S	S	S	S	S	S	S	S	S	S
22	22	22	FR	C	C	C	C	C	C	C	C	C	C	C
23	23	23	FR	A	A	A	A	A	A	A	A	A	A	A
24	24	24	FR	A	A	A	A	A	A	A	A	A	A	A
25	25	25	FR	S	S	S	S	S	S	S	S	S	S	S
26	26	26	CDR-H1	G	G	G	G	G	G	G	G	G	G	G
27	27	27	CDR-H1	F	F	F	F	F	F	F	F	F	F	F
28	28	28	CDR-H1	T	T	T	T	T	T	T	T	T	T	T
29	29	29	CDR-H1	F	F	F	F	F	F	F	F	F	F	F
30	30	30	CDR-H1	S	S	S	S	S	S	S	S	S	S	S
31	31	31	CDR-H1	S	S	S	S	S	S	S	S	S	S	S
32	32	32	CDR-H1	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
33	33	33	CDR-H1	T	E	G	T	T	T	T	T	T	T	T
34	34	34	CDR-H1	M	M	M	M	M	M	M	M	M	M	M

55

**Table 7**  
**Humanized 9D5 Vh Regions**

5	Chothia residue #	Kabat residue #	Linear residue #	FR of CDR	Marine 9D5 VH (SEQ ID NO: 1)	Hu VH Acceptor #1 Ace#BAC03114 (SEQ ID NO: 3)	Hu VH Acceptor #2 Ace#AA82494 (SEQ ID NO: 4)	Hu9D5 VHv1 (SEQ ID NO: 5)	Hu9D5 VHv2 (SEQ ID NO: 6)	Hu9D5 VHv2b (SEQ ID NO: 7)	Hu9D5 VHv3 (SEQ ID NO: 8)	Hu9D5 VHv3b (SEQ ID NO: 9)	Hu9D5 VHv4 (SEQ ID NO: 10)	Hu9D5 VHv4b (SEQ ID NO: 11)	Hu9D5 VHv5 (SEQ ID NO: 12)
10	35	35	35	CDR-H1	S	N	S	S	S	S	S	S	S	S	S
15	36	36	36	Fr2	W	W	W	W	W	W	W	W	W	W	W
20	37	37	37	Fr2	V	V	V	V	V	V	V	V	V	V	V
25	38	38	38	Fr2	R	R	R	R	R	R	R	R	R	R	R
30	39	39	39	Fr2	Q	Q	Q	Q	Q	Q	Q	Q	Q	Q	Q
35	40	40	40	Fr2	T	A	T	A	A	T	A	T	A	A	T
40	41	41	41	Fr2	F	P	P	P	P	P	F	P	P	P	P
45	42	42	42	Fr2	E	G	D	G	G	E	G	G	G	G	E
50	43	43	43	Fr2	K	K	K	K	K	K	K	K	K	K	K
55	44	44	44	Fr2	R	G	R	G	G	R	G	R	G	R	R
60	45	45	45	Fr2	L	L	L	L	L	L	L	L	L	L	L
65	46	46	46	Fr2	E	E	E	E	E	E	E	E	E	E	E
70	47	47	47	Fr2	L	W	W	L	L	L	W	L	W	L	L
75	48	48	48	Fr2	V	V	V	V	V	V	V	V	V	V	V
80	49	49	49	Fr2	A	S	A	A	A	A	S	A	S	A	A
85	50	50	50	CDR-H2	E	Y	I	E	E	E	E	E	E	E	E
90	51	51	51	CDR-H2	I	I	I	I	I	I	I	I	I	I	I
95	52	52	52	CDR-H2	S	I	S	S	S	S	S	S	S	S	S
100	52A	52A	53	CDR-H2	N	S	S	N	N	N	N	N	N	N	N
105	53	53	54	CDR-H2	S	S	G	S	S	S	S	S	S	S	S
110	54	54	55	CDR-H2	G	G	G	G	G	G	G	G	G	G	G
115	55	55	56	CDR-H2	D	S	S	D	D	D	D	D	D	D	D
120	56	56	57	CDR-H2	T	T	Y	T	T	T	T	T	T	T	T
125	57	57	58	CDR-H2	T	I	T	T	T	T	T	T	T	T	T
130	58	58	59	CDR-H2	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
135	59	59	60	CDR-H2	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
140	60	60	61	CDR-H2	P	A	P	P	P	P	P	P	P	P	P
145	61	61	62	CDR-H2	D	D	D	D	D	D	D	D	D	D	D

55

**Table 7**  
**Humanized 9D5 Vh Regions**

5  
10  
15  
20  
25  
30  
35  
40  
45  
50  
55

Chothia residue #	Kabat residue #	Linear residue #	FR of CDR	Marine 9D5 VH (SEQ ID NO: 1)	Hu VH Acceptor Fr Ace#BAC03114 (SEQ ID NO: 3)	Hu VH Acceptor Fr Ace#AA82494 (SEQ ID NO: 4)	Hu9D5 VHv1 (SEQ ID NO: 5)	Hu9D5 VHv2 (SEQ ID NO: 6)	Hu9D5 VHv2b (SEQ ID NO: 7)	Hu9D5 VHv3 (SEQ ID NO: 8)	Hu9D5 VHv3b (SEQ ID NO: 9)	Hu9D5 VHv4 (SEQ ID NO: 10)	Hu9D5 VHv4b (SEQ ID NO: 11)	Hu9D5 VHv5 (SEQ ID NO: 12)
62	62	63	CDR-H2	T	S	S	T	T	T	T	T	T	T	T
63	63	64	CDR-H2	V	V	V	V	V	V	V	V	V	V	V
64	64	65	CDR-H2	R	R	R	R	R	R	R	R	R	R	R
65	65	66	CDR-H2	G	G	G	G	G	G	G	G	G	G	G
66	66	67	Fr3	R	R	R	R	R	R	R	R	R	R	R
67	67	68	Fr3	F	F	F	F	F	F	F	F	F	F	F
68	68	69	Fr3	T	T	T	T	T	T	T	T	T	T	T
69	69	70	Fr3	F	I	I	F	F	I	F	I	I	I	I
70	70	71	Fr3	S	S	S	S	S	S	S	S	S	S	S
71	71	72	Fr3	R	R	R	R	R	R	R	R	R	R	R
72	72	73	Fr3	D	D	D	D	D	D	D	D	D	D	D
73	73	74	Fr3	N	N	N	N	N	N	N	N	N	N	N
74	74	75	Fr3	A	A	A	A	A	A	A	A	A	A	A
75	75	76	Fr3	K	K	K	K	K	K	K	K	K	K	K
76	76	77	Fr3	N	N	N	N	N	N	N	N	N	N	N
77	77	78	Fr3	T	S	T	S	S	T	S	T	S	T	T
78	78	79	Fr3	L	L	L	L	L	L	L	L	L	L	L
79	79	80	Fr3	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
80	80	81	Fr3	L	L	L	L	L	L	L	L	L	L	L
81	81	82	Fr3	Q	Q	Q	Q	Q	Q	Q	Q	Q	Q	Q
82	82	83	Fr3	M	M	M	S	S	M	S	M	S	M	M
82A	82A	84	Fr3	S	N	S	N	N	S	N	S	N	S	N
82B	82B	85	Fr3	S	S	S	S	L	S	L	S	L	S	L
82C	82C	86	Fr3	L	L	L	L	L	L	L	L	L	L	L
83	83	87	Fr3	K	R	K	K	R	K	R	K	R	K	R
84	84	88	Fr3	S	A	S	A	A	S	A	S	A	S	A
85	85	89	Fr3	E	E	E	E	E	E	E	E	E	E	E
86	86	90	Fr3	D	D	D	D	D	D	D	D	D	D	D

**Table 7**  
**Humanized 9D5 Vh Regions**

Chothia residue #	Kabat residue #	Linear residue #	FR of CDR	Marine 9D5 VH (SEQ ID NO: 1)	Hu VH Acceptor #1 Ace#EAC03114 (SEQ ID NO: 3)	Hu VH Acceptor #r Ace#AA82494 (SEQ ID NO: 4)	Hu9D5 VHv1 (SEQ ID NO: 5)	Hu9D5 VHv2 (SEQ ID NO: 6)	Hu9D5 VHv2b (SEQ ID NO: 7)	Hu9D5 VHv3 (SEQ ID NO: 8)	Hu9D5 VHv3b (SEQ ID NO: 9)	Hu9D5 VHv4 (SEQ ID NO: 10)	Hu9D5 VHv4b (SEQ ID NO: 11)	Hu9D5 VHv5 (SEQ ID NO: 12)
87	87	91	Fr3	T	T	T	T	T	T	T	T	T	T	T
88	88	92	Fr3	A	A	A	A	A	A	A	A	A	A	A
89	89	93	Fr3	M	V	M	V	V	M	V	M	V	M	V
90	90	94	Fr3	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
91	91	95	Fr3	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
92	92	96	Fr3	C	C	C	C	C	C	C	C	C	C	C
93	93	97	Fr3	A	A	A	A	A	A	A	A	A	A	A
94	94	98	Fr3	R	R	R	R	R	R	R	R	R	R	R
95	95	99	CDR-H3	H	G	L	H	H	H	H	H	H	H	H
96	96	100	CDR-H3	Y	G	Y	Y	Y	Y	Y	Y	Y	Y	Y
97	97	101	CDR-H3	Y	Q	Y	Y	Y	Y	Y	Y	Y	Y	Y
98	98	102	CDR-H3	Y	G	G	Y	Y	Y	Y	Y	Y	Y	Y
99	99	103	CDR-H3	G	S	Y	G	G	G	G	G	G	G	G
100	100	104	CDR-H3	G	R	R	G	G	G	G	G	G	G	G
100A	100A	105	CDR-H3	G	Y	Y	G	G	G	G	G	G	G	G
100B	100B	106	CDR-H3	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
100C	100C	107	CDR-H3	G	Y	F	G	G	G	G	G	G	G	G
100D	100D	108	CDR-H3	G	Y	-	G	G	G	G	G	G	G	G
100E	100E	109	CDR-H3	W	Y	-	W	W	W	W	W	W	W	W
100F	100F	110	CDR-H3	F	Q	-	F	F	F	F	F	F	F	F
100G	100G	111	CDR-H3	F	M	-	F	F	F	F	F	F	F	F
101	101	112	CDR-H3	D	D	D	D	D	D	D	D	D	D	D
102	102	113	CDR-H3	V	V	Y	V	V	V	V	V	V	V	V
103	103	114	Fr4	W	W	W	W	W	W	W	W	W	W	W
104	104	115	Fr4	G	G	G	G	G	G	G	G	G	G	G
105	105	116	Fr4	T	Q	Q	Q	Q	Q	Q	Q	Q	Q	Q
106	106	117	Fr4	G	G	G	G	G	G	G	G	G	G	G
107	107	118	Fr4	T	T	T	T	T	T	T	T	T	T	T

55

**Table 7**  
**Humanized 9D5 Vh Regions**

Chothia residue #	Kabat residue #	Linear residue #	FR or CDR	Murine 9D5 VH (SEQ ID NO: 1)	Hu VH Acceptor Fr Acc#BAC02114 (SEQ ID NO: 3)	Hu VH Acceptor Fr Acc#AAV82404 (SEQ ID NO: 4)	Hu9D5 VHV1 (SEQ ID NO: 5)	Hu9D5 VHV2 (SEQ ID NO: 6)	Hu9D5 VHV2b (SEQ ID NO: 7)	Hu9D5 VHV3 (SEQ ID NO: 8)	Hu9D5 VHV3b (SEQ ID NO: 9)	Hu9D5 VHV4 (SEQ ID NO: 10)	Hu9D5 VHV4b (SEQ ID NO: 11)	Hu9D5 VHV5 (SEQ ID NO: 12)
108	108	119	Fr4	T	V	M	T	T	L	T	L	T	L	T
109	109	120	Fr4	V	V	V	V	V	V	V	V	V	V	V
110	110	121	Fr4	T	V	V	T	T	T	T	V	T	T	T
111	111	122	Fr4	V	V	V	V	V	V	V	V	V	V	V
112	112	123	Fr4	S	S	S	S	S	S	S	S	S	S	S
113	113	124	Fr4	S	S	S	S	S	S	S	S	S	S	S

**Table 8**  
**Humanized 9D5 Vh Regions**

Chothia residue #	Kabat residue #	Linear residue #	FR or CDR	Murine 9D5 VL (SEQ ID NO: 16)	Hu VL Acceptor Fr Acc#ABC66852 (SEQ ID NO: 18)	Hu9D5 VLV1 (SEQ ID NO: 19)	Hu9D5 VLV2 (SEQ ID NO: 20)	Hu9D5 VLV3 (SEQ ID NO: 21)	Hu9D5 VLV4 (SEQ ID NO: 22)	Hu9D5 VLV5 (SEQ ID NO: 23)
1	1	1	Fr1	D	D	D	D	D	D	D
2	2	2	Fr1	I	V	I	I	I	I	I
3	3	3	Fr1	V	V	V	V	V	V	V
4	4	4	Fr1	M	M	M	M	M	M	M
5	5	5	Fr1	T	V	T	T	T	T	T
6	6	6	Fr1	Q	Q	Q	Q	Q	Q	Q
7	7	7	Fr1	A	S	S	S	S	S	S
8	8	8	Fr1	A	P	P	P	P	A	A
9	9	9	Fr1	F	L	L	L	L	P	P
10	10	10	Fr1	S	S	S	S	S	S	S
11	11	11	Fr1	V	L	L	L	L	L	L
12	12	12	Fr1	F	P	P	P	P	P	P

**Table 8**  
**Humanized 9D5 Vk Regions**

	Chothia residue #	Kabat residue #	Linear residue #	FR of CDR	Mouse 9D5 VL (SEQ ID NO: 16)	Hum VL Acceptor Fr Acc#ABC66952 (SEQ ID NO: 18)	Hum9D5 VLv1 (SEQ ID NO: 19)	Hum9D5 VLv2 (SEQ ID NO: 20)	Hum9D5 VLv3 (SEQ ID NO: 21)	Hum9D5 VLv4 (SEQ ID NO: 22)	Hum9D5 VLv5 (SEQ ID NO: 23)
5											
10											
15											
20											
25											
30											
35											
40											
45											
50											
55											
	13	13	13	Fr1	V	V	V	V	V	V	
	14	14	14	Fr1	T	T	T	T	T	T	
	15	15	15	Fr1	F	F	F	F	F	F	
	16	16	16	Fr1	G	G	G	G	G	G	
	17	17	17	Fr1	E	E	E	E	E	E	
	18	18	18	Fr1	S	P	P	P	P	S	
	19	19	19	Fr1	V	A	A	A	V	V	
	20	20	20	Fr1	S	S	S	S	S	S	
	21	21	21	Fr1	I	I	I	I	I	I	
	22	22	22	Fr1	S	S	S	S	S	S	
	23	23	23	Fr1	C	C	C	C	C	C	
	24	24	24	CDR-L1	R	R	R	R	R	R	
	25	25	25	CDR-L1	S	S	S	S	S	S	
	26	26	26	CDR-L1	S	N	S	S	S	S	
	27	27	27	CDR-L1	K	Q	K	K	K	K	
	27A	27A	28	CDR-L1	S	S	S	S	S	S	
	27B	27B	29	CDR-L1	L	L	L	L	L	L	
	27C	27C	30	CDR-L1	L	L	L	L	L	L	
	27D	27D	31	CDR-L1	H	Y	H	H	H	H	
	27E	27E	32	CDR-L1	S	S	S	S	S	S	
	28	28	33	CDR-L1	N	N	N	N	N	N	
	29	29	34	CDR-L1	G	G	G	G	G	G	
	30	30	35	CDR-L1	N	Y	N	N	N	N	
	31	31	36	CDR-L1	T	N	T	T	T	T	
	32	32	37	CDR-L1	Y	Y	Y	Y	Y	Y	
	33	33	38	CDR-L1	L	L	L	L	L	L	
	34	34	39	CDR-L1	Y	D	Y	Y	Y	Y	
	35	35	40	Fr2	W	W	W	W	W	W	
	36	36	41	Fr2	F	Y	F	Y	F	F	
	37	37	42	Fr2	L	L	L	L	L	L	

**Table 8**  
**Humanized 9D5 Vk Regions**

5  
10  
15  
20  
25  
30  
35  
40  
45  
50  
55

Chothia residue #	Kabat residue #	Linear residue #	FR of CDR	Mouse 9D5 VL (SEQ ID NO: 16)	Hum VL Acceptor Fr Acc#ABC66952 (SEQ ID NO: 18)	Hum9D5 VLv1 (SEQ ID NO: 19)	Hum9D5 VLv2 (SEQ ID NO: 20)	Hum9D5 VLv3 (SEQ ID NO: 21)	Hum9D5 VLv4 (SEQ ID NO: 22)	Hum9D5 VLv5 (SEQ ID NO: 23)
38	38	43	Fr2	Q	Q	Q	Q	Q	Q	Q
39	39	44	Fr2	R	R	R	R	R	R	R
40	40	45	Fr2	F	F	F	F	F	F	F
41	41	46	Fr2	G	G	G	G	G	G	G
42	42	47	Fr2	Q	Q	Q	Q	Q	Q	Q
43	43	48	Fr2	S	S	S	S	S	S	S
44	44	49	Fr2	F	F	F	F	F	F	F
45	45	50	Fr2	Q	Q	Q	Q	Q	Q	Q
46	46	51	Fr2	L	L	L	L	L	L	L
47	47	52	Fr2	L	L	L	L	L	L	L
48	48	53	Fr2	I	I	I	I	I	I	I
49	49	54	Fr2	Y	Y	Y	Y	Y	Y	Y
50	50	55	CDR-L2	R	S	R	R	R	R	R
51	51	56	CDR-L2	V	G	V	V	V	V	V
52	52	57	CDR-L2	S	S	S	S	S	S	S
53	53	58	CDR-L2	N	N	N	N	N	N	N
54	54	59	CDR-L2	L	R	L	L	L	L	L
55	55	60	CDR-L2	A	A	A	A	A	A	A
56	56	61	CDR-L2	S	S	S	S	S	S	S
57	57	62	Fr3	G	G	G	G	G	G	G
58	58	63	Fr3	V	V	V	V	V	V	V
59	59	64	Fr3	F	P	F	F	F	F	F
60	60	65	Fr3	D	D	D	D	S	S	S
61	61	66	Fr3	R	R	R	R	R	R	R
62	62	67	Fr3	F	F	F	F	F	F	F
63	63	68	Fr3	S	S	S	S	S	S	S
64	64	69	Fr3	G	G	G	G	G	G	G
65	65	70	Fr3	S	S	S	S	S	S	S
66	66	71	Fr3	G	G	G	G	G	G	G
67	67	72	Fr3	S	S	S	S	S	S	S

**Table 8**  
**Humanized 9D5 Vk Regions**

	Chothia residue #	Kabat residue #	Linear residue #	FR of CDR	Mouse 9D5 VL (SEQ ID NO: 16)	Hum VL Acceptor Fr Acc#ABC66952 (SEQ ID NO: 18)	Hum9D5 VLv1 (SEQ ID NO: 19)	Hum9D5 VLv2 (SEQ ID NO: 20)	Hum9D5 VLv3 (SEQ ID NO: 21)	Hum9D5 VLv4 (SEQ ID NO: 22)	Hum9D5 VLv5 (SEQ ID NO: 23)
5											
10											
15											
20											
25											
30											
35											
40											
45											
50											
55											
	68	68	73	Fc3	G	G	G	G	G	G	G
	69	69	74	Fc3	T	T	T	T	T	T	T
	70	70	75	Fc3	A	D	D	D	D	A	A
	71	71	76	Fc3	F	F	F	F	F	F	F
	72	72	77	Fc3	T	T	T	T	T	T	T
	73	73	78	Fc3	L	L	L	L	L	L	L
	74	74	79	Fc3	R	K	K	K	K	R	R
	75	75	80	Fc3	I	I	I	I	I	I	I
	76	76	81	Fc3	S	S	S	S	S	S	S
	77	77	82	Fc3	R	R	R	R	R	R	R
	78	78	83	Fc3	V	V	V	V	V	V	V
	79	79	84	Fc3	E	E	E	E	E	E	E
	80	80	85	Fc3	A	A	A	A	A	A	A
	81	81	86	Fc3	E	E	E	E	E	E	E
	82	82	87	Fc3	D	D	D	D	D	D	D
	83	83	88	Fc3	V	V	V	V	V	V	V
	84	84	89	Fc3	G	G	G	G	G	G	G
	85	85	90	Fc3	V	V	V	V	V	V	V
	86	86	91	Fc3	Y	Y	Y	Y	Y	Y	Y
	87	87	92	Fc3	Y	Y	Y	Y	Y	Y	Y
	88	88	93	Fc3	C	C	C	C	C	C	C
	89	89	94	CDR-L3	M	M	M	M	M	M	M
	90	90	95	CDR-L3	Q	Q	Q	Q	Q	Q	Q
	91	91	96	CDR-L3	H	A	H	H	H	H	H
	92	92	97	CDR-L3	L	L	L	L	L	L	L
	93	93	98	CDR-L3	E	Q	E	E	E	E	E
	94	94	99	CDR-L3	Y	S	Y	Y	Y	Y	Y
	95	95	100	CDR-L3	P	P	P	P	P	P	P
	96	96	101	CDR-L3	L	Y	L	L	L	L	L
	97	97	102	CDR-L3	T	T	T	T	T	T	T

**Table 8**  
**Humanized 9D5 VK Regions**

5	Chothia residue #	Kabat residue #	Linear residue #	FR or CDR	Mouse 9D5 VL (SEQ ID NO: 16)	Hu VL Acceptor Fr Acc#ABC66952 (SEQ ID NO: 18)	Hu9D5 VLv1 (SEQ ID NO: 19)	Hu9D5 VLv2 (SEQ ID NO: 20)	Hu9D5 VLv3 (SEQ ID NO: 21)	Hu9D5 VLv4 (SEQ ID NO: 22)	Hu9D5 VLv5 (SEQ ID NO: 23)
10	98	98	103	F <sub>H</sub>	F	F	F	F	F	F	F
15	99	99	104	F <sub>H</sub>	G	G	G	G	G	G	G
	100	100	105	F <sub>H</sub>	A	Q	Q	Q	Q	Q	Q
	101	101	106	F <sub>H</sub>	G	G	G	G	G	G	G
	102	102	107	F <sub>H</sub>	T	T	T	T	T	T	T
20	103	103	108	F <sub>H</sub>	K	K	K	K	K	K	K
	104	104	109	F <sub>H</sub>	L	L	L	L	L	L	L
	105	105	110	F <sub>H</sub>	E	E	E	E	E	E	E
25	106	106	111	F <sub>H</sub>	L	I	I	I	I	I	I
	107	107	112	F <sub>H</sub>	K	K	K	K	K	K	K

**Table 9**

<b>V<sub>H</sub>, V<sub>L</sub> Backmutations and Other Mutations</b>		
<b>V<sub>H</sub> or V<sub>L</sub> Variant</b>	<b>V<sub>H</sub> or V<sub>L</sub> Exon Acceptor Sequence</b>	<b>Donor Framework Residues</b>
Hu9D5VHv1 (SEQ ID NO: 5)	NCBI accession codes BAC02114 and AAX8294 (SEQ ID NOS: 3 and 4, respectively)	H47, H69, H82
35 Hu9D5VHv2 (SEQ ID NO: 6)	NCBI accession codes BAC02114 and AAX8294 (SEQ ID NOS: 3 and 4, respectively)	H47, H69, H82, H82b
Hu9D5VHv2b (SEQ ID NO: 7)	NCBI accession codes BAC02114 and AAX8294 (SEQ ID NOS: 3 and 4, respectively)	H42, H47, H108
40 Hu9D5VHv3 (SEQ ID NO: 8)	NCBI accession codes BAC02114 and AAX8294 (SEQ ID NOS: 3 and 4, respectively)	H69, H82, H82b
Hu9D5VHv3b (SEQ ID NO: 9)	NCBI accession codes BAC02114 and AAX8294 (SEQ ID NOS: 3 and 4, respectively)	H47, H108
45 Hu9D5VHv4 (SEQ ID NO: 10)	NCBI accession codes BAC02114 and AAX8294 (SEQ ID NOS: 3 and 4, respectively)	H82, H82b
Hu9D5VHv4b (SEQ ID NO: 11)	NCBI accession codes BAC02114 and AAX8294 (SEQ ID NOS: 3 and 4, respectively)	H47, H108
50 Hu9D5VHv5 (SEQ ID NO: 12)	NCBI accession codes BAC02114 and AAX8294 (SEQ ID NOS: 3 and 4, respectively)	H42, H47, H82b
Hu9D5VLv1 (SEQ ID NO: 19)	NCBI accession code ABC66952 (SEQ ID NO: 18)	L36
55 Hu9D5VLv2 (SEQ ID NO: 20)	NCBI accession code ABC66952 (SEQ ID NO: 18)	None

(continued)

<b>V<sub>H</sub>, V<sub>L</sub> Backmutations and Other Mutations</b>		
V <sub>H</sub> or V <sub>L</sub> Variant	V <sub>H</sub> or V <sub>L</sub> Exon Acceptor Sequence	Donor Framework Residues
Hu9D5VLv3 (SEQ ID NO: 21)	NCBI accession code ABC66952 (SEQ ID NO: 18)	L60
Hu9D5VLv4 (SEQ ID NO: 22)	NCBI accession code ABC66952 (SEQ ID NO: 18)	L8, L9, L19, L36, L39, L60, L70, L74
Hu9D5VLv5 (SEQ ID NO: 23)	NCBI accession code ABC66952 (SEQ ID NO: 18)	L8, L9, L18, L19, L36, L39, L60, L70, L74

**Table 10**  
**Kabat Numbering of Framework Residues**  
**for Backmutations and Other Mutations in V<sub>H</sub> Regions of Humanized 9D5 Antibodies**

Residue	BAC02114 Heavy Chain	AAX82494 Heavy Chain	Mouse 9D5	Hu9D5 VHv1	Hu9D5 VHv2	Hu9D5 VHv2b	Hu9D5 VHv3	Hu9D5 VHv3b	Hu9D5 VHv4	Hu9D5 VHv4b	Hu9D5 VHv5
H42	G	D	E	G	G	E	G	G	G	G	E
H47	W	W	L	L	L	L	W	L	W	L	L
H69	I	I	F	F	F	I	F	I	I	I	I
H82	M	M	M	S	S	M	S	M	S	M	M
H82b	S	S	S	S	L	S	L	S	L	S	L
H108	T	T	T	T	T	L	T	L	T	L	T

**Table 11**

<b>Kabat Numbering of Framework Residues for Backmutations and Other Mutations in V<sub>K</sub> Regions of Humanized 9D5 Antibodies</b>							
Residue	ABC66952 Light Chain	Mouse 9D5	Hu9D5 VLv1	Hu9D5 VLv2	Hu9D5 VLv3	Hu9D5 VLv4	Hu9D5 VLv5
L8	P	A	P	P	P	A	A
L9	L	P	L	L	L	P	P
L18	P	s	P	P	P	P	s
L19	A	V	A	A	A	V	V
L36	Y	F	F	Y	Y	F	F
L39	K	R	K	K	K	R	R
L60	D	D	D	D	S	S	S
L70	D	A	D	D	D	A	A
L74	K	R	K	K	K	R	R

**[0227]** An alignment of the murine 9D5 V<sub>H</sub> sequence (SEQ ID NO: 1) with the mouse model sequences (1SEQ\_H and 1MQK\_H; SEQ ID NOS: 2 and 62, respectively), the human acceptor sequences (BAC02114 and AAX82494; SEQ ID NOS: 3 and 4, respectively), and the Hu9D5VHv1, Hu9D5VHv2, Hu9D5VHv2b, Hu9D5VHv3, Hu9D5VHv3b,

Hu9D5VHv4, Hu9D5VHv4b, and Hu9D5VHv5 sequences (SEQ ID NOS: 5-12, respectively), is shown in Figure 1A. The CDRs as defined by Kabat are enclosed in boxes, except that the first enclosed box is a composite of the Chothia CDR-H1 and the Kabat CDR-H1, with the Kabat CDR-H1 underlined and bolded. Positions at which canonical, vernier, or interface residues differ between mouse and human acceptor sequences are candidates for substitution. Examples of vernier/CDR foundation residues include residues 2, 49, 69, 71, 75, 80, and 94 by Kabat numbering in Table 7. Examples of canonical/CDR interacting residues include residues 24, 48, and 73 by Kabat numbering in Table 7. Examples of interface/packing (VH+VL) residues include residues 37, 39, 44, 47, 91, 93, and 103 by Kabat numbering in Table 7.

**[0228]** An alignment of the murine 9D5 Vk sequence (SEQ ID NO: 16) with the mouse model sequence (1MJU L; SEQ ID NO: 17), the human acceptor sequence (ABC66952; SEQ ID NO: 18), and the Hu9D5VLv1, Hu9D5VLv2, Hu9D5VLv3, Hu9D5VLv4, and Hu9D5VLv5 sequences (SEQ ID NOS: 19-23, respectively), is shown in Figure 1B. The CDRs as defined by Kabat are enclosed in boxes, except that the first enclosed box is a composite of the Chothia CDR-H1 and the Kabat CDR-H1, with the Kabat CDR-H1 underlined and bolded. Positions at which canonical, vernier, or interface residues differ between mouse and human acceptor sequences are candidates for substitution. Examples of vernier/CDR foundation residues include residues 4, 35, 46, 49, 66, 68, and 69 by Kabat numbering in Table 8. Examples of canonical/CDR interacting residues include residues 2, 48, 64, and 71 by Kabat numbering in Table 8. Examples of interface/packing (VH+VL) residues include residues 36, 38, 44, 47, 87, and 98 by Kabat numbering in Table 8.

**[0229]** The rationales for selection of the positions indicated in Tables 9 and 11 in the light chain variable region as candidates for substitution are as follows.

**[0230]** P8A: The model shows a kink in the loop at this position, so a backmutation to A was tried to alleviate this.

**[0231]** L9P: The model shows a kink in the loop at this position, so a backmutation to P was tried to alleviate this.

**[0232]** P18S: P at this position is less frequent. Mutation to S was tried to alleviate loop distortion.

**[0233]** A19V: A and V at this position are equally frequent. Mutation to V was tried to alleviate loop distortion.

**[0234]** Y36F: This is an interface residue, and typically Y or F. Y has an extra hydroxyl group that could potentially affect LC + HC packing. Polarity of the Y36 could potentially interfere with placement of heavy chain CRD-H3 residue H95. A homology model shows that Y at this position will form a de-novo H-bond with F100(g) in H3 that may result in mobility restriction of H3. Both F and Y were used in separate versions.

**[0235]** K39R: R forms H-bonds with adjoining residues in this loop as compared to K. To rule out any effects on loop stability, mutation to R was tried.

**[0236]** D60S: The presence of D at this residue shows high exposure for proteolysis. In some versions it was replaced with S, a residue most frequent in human germ line at this position. This is predicted to enhance stability.

**[0237]** D70A: D has proteolysis potential, so mutation to A was tried.

**[0238]** K74R: R appears to stabilize the loop compared to K at this position, so backmutation to R was tried.

**[0239]** The rationales for selection of the positions indicated in Tables 9 and 10 in the heavy chain variable region as candidates for substitution are as follows.

**[0240]** G42E: E makes ionic interactions with R at position 44. To rule out any effect that these interactions might have, backmutation to E was tested in some versions.

**[0241]** W47L: This is an interface residue, typically W. In murine 9D5 heavy chain, it is L, whereas in the human acceptor framework there is W at this position. Although L and W are both non-polar, the phenolic ring in W could potentially impact light chain: heavy chain packing. W and L were included in separate versions.

**[0242]** I69F: This is a vernier residue, part of the CDR-H2 foundation. According to the homology model, the aromatic ring of F in murine sequences makes a pie stack with the aromatic ring of CRD-H2 residue Y59. An I residue at this position disturbs that stacking. I and F were included in separate versions.

**[0243]** M82S: M at this position is very rare in human frameworks. More common are A or N. Changing the residue to the more common S might reduce immunogenicity that could be associated with the rare M at this position. Based upon model observations, there is a possibility that M interacts with Leu80, which is a vernier zone residue. M and S were included in separate versions.

**[0244]** S82(b)L: S is present at this position in both the murine and human framework sequences, but S at this position is less frequent. Judging from the position where this residue sits in the model, it possibly could make contact with the antigen. Residues in framework region 3 of the heavy chain are known to occasionally contribute towards binding. S and L were included in separate versions.

**[0245]** T108L: L is the most frequent residue in human frameworks at this position, so L was tried in some versions.

**[0246]** Because two human acceptor frameworks (BAC02114 and AAX82494.1 (SEQ ID NOS: 3 and 4)) were used for humanization of the 9D5 heavy chain variable region, there are certain framework positions that have different amino acids in the two acceptor frameworks. Humanized versions of the 9D5 heavy chain variable region can include either amino acid at these residues. Examples of residues in which these two acceptors differ include positions H19 (R or K), H40 (A or T), H44 (G or R), H49 (S or A), H77 (S or T), H82a (N or S), H83 (R or K), H84 (A or S), and H89 (V or M) by Kabat numbering. The rationales for choosing one amino acid or the other at these positions are as follows.

**[0247]** H19 (R or K): R and K are present in the two frameworks considered, so each residue was tried in separate

versions.

**[0248]** H40 (A or T): A and T are present in the two frameworks considered, so each residue was tried in separate versions.

**[0249]** H44 (G or R): G and R are present in the two frameworks considered, so each residue was tried in separate versions.

**[0250]** H49 (S or A): This is a vernier residue that packs under CDR-H2. In the murine sequence, it is A. The slight size difference and hydrophilic nature of S could potentially perturb the CDR foundation. S and A were included in separate versions.

**[0251]** H77 (S or T): S and T are present in the two human acceptor frameworks considered, so each was tried in separate versions.

**[0252]** H82a (N or S): N at this position is much less frequent than S. Moreover, S appears to be contributing to antigen binding. Mutation to S was tried in some versions.

**[0253]** H83 (R or K): K in framework region 3 is close to the CDR binding surface area. Replacing it with R, which is bulkier than K, might obstruct placement of antigen. R is most frequent and K is third frequent in human frameworks at this position. R and K were included in separate versions.

**[0254]** H84 (A or S): S and A are present in the two frameworks considered, so each residue was tried in separate versions.

**[0255]** H89 (V or M): V and M are present in the two frameworks considered, so each was tried in separate versions.

**[0256]** The five humanized light chain variable region variants and five humanized heavy chain variable region variants are as follows:

**[0257]** Hu9D5VL version 1 (Y36F substitution in lowercase):

DIVMTQSP<sub>1</sub>SLPVT<sub>1</sub>PGEPASISCRSSK<sub>1</sub>SLLSHNSGNTYLYWFLQKPGQSPQLLIYR<sub>1</sub>VSNLAS  
GVPDRFSGSGSGTDFTLKISRVEAEDVGVYYCMQHLEYPLTFGQGTKLEIK (SEQ ID NO:  
19)

**[0258]** Hu9D5VL version 2 (no substitutions):

DIVMTQSP<sub>1</sub>SLPVT<sub>1</sub>PGEPASISCRSSK<sub>1</sub>SLLSHNSGNTYLYWYFLQKPGQSPQLLIYR<sub>1</sub>VSNLAS  
GVPDRFSGSGSGTDFTLKISRVEAEDVGVYYCMQHLEYPLTFGQGTKLEIK (SEQ ID NO:  
20)

**[0259]** Hu9D5VL version 3 (D60S substitution in lowercase):

DIVMTQSP<sub>1</sub>SLPVT<sub>1</sub>PGEPASISCRSSK<sub>1</sub>SLLSHNSGNTYLYWYFLQKPGQSPQLLIYR<sub>1</sub>VSNLAS  
GVP<sub>1</sub>RFSGSGSGTDFTLKISRVEAEDVGVYYCMQHLEYPLTFGQGTKLEIK (SEQ ID NO:  
21)

**[0260]** Hu9D5VL version 4 (P8A, L9P, A19V, Y36F, K39R, D60S, D70A, and K74R substitutions in lowercase):

DIVMTQ<sub>1</sub>S<sub>1</sub>apSLPVT<sub>1</sub>PGEP<sub>1</sub>SISCRSSK<sub>1</sub>SLLSHNSGNTYLYWFLQ<sub>1</sub>PGQSPQLLIYR<sub>1</sub>VSNLASG  
VP<sub>1</sub>RFSGSGSGT<sub>1</sub>aFTL<sub>1</sub>ISRVEAEDVGVYYCMQHLEYPLTFGQGTKLEIK (SEQ ID NO: 22)

**[0261]** Hu9D5VL version 5 (P8A, L9P, P18S, A19V, Y36F, K39R, D60S, D70A, and K74R substitutions in lowercase):

DIVMTQ<sub>1</sub>S<sub>1</sub>apSLPVT<sub>1</sub>PGE<sub>1</sub>svSISCRSSK<sub>1</sub>SLLSHNSGNTYLYWFLQ<sub>1</sub>PGQSPQLLIYR<sub>1</sub>VSNLASGV  
P<sub>1</sub>RFSGSGSGT<sub>1</sub>aFTL<sub>1</sub>ISRVEAEDVGVYYCMQHLEYPLTFGQGTKLEIK (SEQ ID NO: 23)

**[0262]** Hu9D5VH version 1 (W47L, I69F, and M82S substitutions in lowercase):

EVQLVESGGGLVQPGGSLRLSCAASGFTFSSYTMSWVRQAPGKGLEIVAEISNSGDTTY  
 YPDIVKGRFTISRDNKNSLYLQ<sub>s</sub>NSLKAEDTAVYYCARHYYYGGGYGGWFFDVGQ  
 GTITVTVSS (SEQ ID NO: 5)

5

**[0263]** Hu9D5VH version 2 (W47L, I69F, M82S, and S82(b)L substitutions in lowercase):

EVQLVESGGGLVQPGGSLRLSCAASGFTFSSYTMSWVRQAPGKGLEIVAEISNSGDTTY  
 YPDIVKGRFTISRDNKNSLYLQ<sub>s</sub>NILRAEDTAVYYCARHYYYGGGYGGWFFDVGQ  
 GTITVTVSS (SEQ ID NO: 6)

10

**[0264]** Hu9D5VH version 2b (G42E, W47L, and T108L substitutions in lowercase):

EVQLVESGGGLVQPGGSLKLSCAASGFTFSSYTMSWVRQTP<sub>e</sub>KRLEIVAEISNSGDTTY  
 PDIVKGRFTISRDNKNTLYLQ<sub>m</sub>SSLKSEDTAMYVCARHYYYGGGYGGWFFDVGQ  
 GTITVTVSS (SEQ ID NO: 7)

15

**[0265]** Hu9D5VH version 3 (I69F, M82S, and S82(b)L substitutions in lowercase):

EVQLVESGGGLVQPGGSLRLSCAASGFTFSSYTMSWVRQAPGKGLEWVSEISNSGDTTY  
 YPDIVKGRFTISRDNKNSLYLQ<sub>s</sub>NILRAEDTAVYYCARHYYYGGGYGGWFFDVGQ  
 GTITVTVSS (SEQ ID NO: 8)

25

**[0266]** Hu9D5VH version 3b (W47L and T108L substitutions in lowercase):

EVQLVESGGGLVQPGGSLKLSCAASGFTFSSYTMSWVRQTPGKRLEIVAEISNSGDTTY  
 YPDIVKGRFTISRDNKNTLYLQ<sub>m</sub>SSLKSEDTAMYVCARHYYYGGGYGGWFFDVGQ  
 QGTITVTVSS (SEQ ID NO: 9)

30

**[0267]** Hu9D5VH version 4 (M82S and S82(b)L substitutions in lowercase):

EVQLVESGGGLVQPGGSLRLSCAASGFTFSSYTMSWVRQAPGKGLEWVSEISNSGDTTY  
 YPDIVKGRFTISRDNKNSLYLQ<sub>s</sub>NILRAEDTAVYYCARHYYYGGGYGGWFFDVGQ  
 GTITVTVSS (SEQ ID NO: 10)

35

40

**[0268]** Hu9D5VH version 4b (W47L and T108L substitutions in lowercase):

EVQLVESGGGLVQPGGSLKLSCAASGFTFSSYTMSWVRQAPGKRLEIVAEISNSGDTTY  
 YPDIVKGRFTISRDNKNTLYLQ<sub>m</sub>SSLKSEDTAMYVCARHYYYGGGYGGWFFDVGQ  
 QGTITVTVSS (SEQ ID NO: 11)

45

**[0269]** Hu9D5VH version 5 (G42E, W47L, and S82(b)L substitutions in lowercase):

EVQLVESGGGLVQPGGSLRLSCAASGFTFSSYTMSWVRQTP<sub>e</sub>KRLEIVAEISNSGDTTY  
 PDIVKGRFTISRDNKNTLYLQ<sub>m</sub>NILRAEDTAVYYCARHYYYGGGYGGWFFDVGQ  
 TTVTVSS (SEQ ID NO: 12)

50

55

**[0270]** Protein quality analysis along with aggregation potential analyses showed no obvious clusters of aggregation-prone residues present in the 9D5 light chain or heavy chain framework regions or CDRs (Wang et al., mAbs 1(3):254-267, 2009).

**Example 8. Binding Kinetic Analysis of Humanized 9D5 Antibodies**

**[0271]** Binding kinetics of binding of all humanized 9D5 variants, murine 9D5, and chimeric 9D5 to recombinant human TTR F87M/L110M were characterized by Biacore, as shown in Table 12. Anti-human (GE Healthcare) was immobilized on sensor chip C5 (lacking dextran chains) via amine coupling following the instructions provided in the GE Healthcare anti-human kit, and mAbs were captured to a level to ensure a maximum binding of analyte of 30-50 RU. Various concentrations of analyte (recombinant human TTR F87M/L110M) were passed over the captured ligand at 50 ul/min in running buffer (HBS + 0.05% P-20, 1 mg/mL BSA) in 3-fold dilutions. For each concentration, the reaction proceeded for a time frame allowing for the higher analyte concentrations to reach equilibrium during association, as well as at least 10% of signal to decay during dissociation. At least one concentration (not the highest or lowest) was run in duplicate. Concentration ranges of analyte were selected based on preliminary experimentation to span at least 10-fold above  $K_D$  to 10-fold below  $K_D$ .

**[0272]** Table 12 summarizes the Biacore association rate ( $k_a$ ), dissociation rate ( $k_d$ ), and binding affinity constant ( $K_D$ ) of the murine, chimeric, and humanized 9D5 variants for recombinant human TTR F87M/L110M.

**Table 12**

<b>Antigen Binding Affinity of 9D5 Antibodies to Human TTR (F87M/L110M)</b>				
<b>Antibody</b>	<b><math>k_a</math> (1/Ms)</b>	<b><math>k_d</math> (1/s)</b>	<b><math>K_D</math> (M)</b>	<b>Rmax</b>
Murine 9D5	2.72E+04	4.93E-04	1.82E-08	31.55
Chimeric 9D5	2.59E+05	3.50E-04	1.35E-09	47.05
Hu9D5 H2L5	1.15E+05	7.85E-04	6.84E-09	57.6
Hu9D5 H3L1	1.52E+04	6.57E-04	4.32E-09	54.93
Hu9D5 H4L1	2.54E+05	5.33E-04	2.09E-09	43.91

**[0273]** These results indicate that the affinity of murine 9D5 for TTR-F87M/L110M ( $K_D = 1.82E-08M$ ) has been slightly improved in the chimeric 9D5 variant ( $K_D = 1.35E-09M$ ). Furthermore, the fully humanized 9D5 variants all have similar affinities in the low nM range. Hu9D5 H4L1 has the strongest affinity ( $K_D = 2.09E-09$ ) of the humanized 9D5 variants tested.

**Example 9. Design of Humanized 14G8 Antibodies**

**[0274]** The starting point or donor antibody for humanization was the mouse antibody 14G8. The heavy chain variable amino acid sequence of mature m14G8 is provided as SEQ ID NO: 61. The light chain variable amino acid sequence of mature m14G8 is provided as SEQ ID NO: 70. The heavy chain CDR1, CDR2, and CDR3 amino acid sequences are provided as SEQ ID NOS: 67-69, respectively (as defined by Kabat). A composite Chothia-Kabat CDR-H1 is provided as SEQ ID NO: 118. The light chain CDR1, CDR2, and CDR3 amino acid sequences are provided as SEQ ID NOS: 77-79, respectively (as defined by Kabat). A variant version of CDR1 is provided as SEQ ID NO: 80. Kabat numbering is used throughout in this Example.

**[0275]** The variable kappa (Vk) of m14G8 belongs to mouse Kabat subgroup 2, which corresponds to human Kabat subgroup 2. The variable heavy (Vh) of m14G8 belongs to mouse Kabat subgroup 3d, which corresponds to Kabat subgroup 1. See Kabat et al. Sequences of Proteins of Immunological Interest, Fifth Edition. NIH Publication No. 91-3242, 1991. The 16-residue CDR-L1 belongs to canonical class 3, the 7-residue CDR-L2 belongs to canonical class 1, and the 9-residue CDR-L3 belongs to canonical class 1 in Vk. See Martin & Thornton, J. Mol. Biol. 263:800-15, 1996. The 10-residue composite Chothia-Kabat CDR-H1 belongs to canonical class 1, and the 17-residue CDR-H2 belongs to canonical class 1. See Martin & Thornton, J Mol. Biol. 263:800-15, 1996. The CDR-H3 has no canonical classes, but the 15-residue loop probably has a kinked base according to the rules of Shirai et al., FEBS Lett. 455:188-97 (1999).

**[0276]** The residues at the interface between the Vk and Vh domains are the ones commonly found except L47, at which position the principle amino acid is usually W.

**[0277]** A search was made over the protein sequences in the PDB database (Deshpande et al., Nucleic Acids Res. 33: D233-7, 2005) to find structures which would provide a rough structural model of 14G8. The crystal structure of Fab with esterase activity (pdb 1MJU) was used as the model for the Vk structure. It was solved at a resolution of 1.22 Å and contains the same canonical structures for CDR-H1 and CDR-H2, and also contains the same length CDR-H3 with a kinked based. An anti-cytochrome C oxidase antibody 7E2 Fv fragment (pdb code 1MQK\_H) was used for the Vh structure. It has a resolution of 1.28 Å and retains the same canonical structure for the loops as 14G8. BioLuminate software (licensed from Schrodinger Inc.) was used to model a rough structure of 14G8.

[0278] A search of the non-redundant protein sequence database from NCBI allowed selection of suitable human frameworks into which to graft the murine CDRs. For Vh, human Ig heavy chains with NCBI accession codes AAD30410.1 and AAX82494.1 were chosen. These share the canonical forms of 14G8 CDR-H1 and H2, and H3 of AAD30410.1 is 15 residues long with a predicted kinked base. For Vk, two human kappa light chains with NCBI accession codes ABA71374.1 and ABC66952.1 were chosen. They have the same canonical classes for LCDRs.

[0279] Three humanized heavy chain variable region variants and three humanized light chain variable region variants were constructed containing different permutations of substitutions (Hu14G8VHv1-3 (SEQ ID NOS: 64-66, respectively) and Hu14G8VLv1-3 (SEQ ID NOS: 74-76, respectively)) (Tables 13 and 14). The exemplary humanized Vh and Vk designs, with backmutations and other mutations based on selected human frameworks, are shown in Tables 13 and 14, respectively. The gray-shaded areas in the first column in Tables 13 and 14 indicate the CDRs as defined by Chothia, and the gray-shaded areas in the remaining columns in Tables 13 and 14 indicate the CDRs as defined by Kabat. SEQ ID NOS: 64-66 and 74-76 contain backmutations and other mutations as shown in Table 15. The amino acids at positions H1, H3, H47, H105, L8, L9, L19, L26, L36, L60, and L70 in Hu14G8VHv1-3 and Hu14G8VLv1-3 are listed in Tables 16 and 17.

**Table 13**  
**Humanized 14G8 Vh Regions**

Chothia residue #	Kabat residue #	Linear residue #	FR of CDR	Murine 14G8 VH (SEQ ID NO: 61)	Hu VH Acceptor Pt Acc# AAD30410.1 (SEQ ID NO: 63)	Hu VH Acceptor Pt Acc# AAX82494.1 (SEQ ID NO: 64)	Hu14G8 VHv1 (SEQ ID NO: 64)	Hu14G8 VHv2 (SEQ ID NO: 65)	Hu14G8 VHv3 (SEQ ID NO: 66)
1	1	1	FR1	E	Q	Q	E	E	E
2	2	2	FR1	V	V	V	V	V	V
3	3	3	FR1	K	Q	Q	K	Q	Q
4	4	4	FR1	L	L	L	L	L	L
5	5	5	FR1	V	V	Q	V	V	V
6	6	6	FR1	E	Q	E	E	E	E
7	7	7	FR1	S	S	S	S	S	S
8	8	8	FR1	Q	Q	G	Q	Q	Q
9	9	9	FR1	G	G	G	G	G	G
10	10	10	FR1	G	G	G	G	G	G
11	11	11	FR1	L	L	L	L	L	L
12	12	12	FR1	V	V	V	V	V	V

**Table 13**  
**Humanized 14G8 Vh Regions**

	Chelonia residue #	Kabat residue #	Linear residue #	FR or CDR	Murine 14G8 VH (SEQ ID NO: 61)	Hu VH Acceptor Fr Ace# AAD30410.1 (SEQ ID NO: 63)	Hu VH Acceptor Fr Ace# AAK82494.1 (SEQ ID NO: 4)	Hu14G8 VHV1 (SEQ ID NO: 64)	Hu14G8 VHV2 (SEQ ID NO: 65)	Hu14G8 VHV3 (SEQ ID NO: 66)
5	13	13	13	Fr1	Q	Q	K	Q	Q	Q
10	14	14	14	Fr1	P	P	P	P	P	P
15	15	15	15	Fr1	G	G	G	G	G	G
	16	16	16	Fr1	G	G	G	G	G	G
20	17	17	17	Fr1	S	S	S	S	S	S
	18	18	18	Fr1	L	S	L	L	L	L
	19	19	19	Fr1	K	K	K	K	K	K
	20	20	20	Fr1	L	L	L	L	L	L
25	21	21	21	Fr1	S	S	S	S	S	S
	22	22	22	Fr1	C	C	C	C	C	C
	23	23	23	Fr1	A	A	A	A	A	A
	24	24	24	Fr1	A	A	A	A	A	A
30	25	25	25	Fr1	S	S	S	S	S	S
	26	26	26	Fr1	G	G	G	G	G	G
	27	27	27	Fr1	F	F	F	F	F	F
	28	28	28	Fr1	Y	Y	Y	Y	Y	Y
35	29	29	29	Fr1	F	F	F	F	F	F
	30	30	30	Fr1	S	S	S	S	S	S
	31	31	31	CDR-H1	S	S	S	S	S	S
	32	32	32	CDR-H1	Y	Y	Y	Y	Y	Y
40	33	33	33	CDR-H1	T	A	G	T	T	T
	34	34	34	CDR-H1	M	M	M	M	M	M
	35	35	35	CDR-H1	S	S	S	S	S	S
45	36	36	36	Fr2	W	W	W	W	W	W
	37	37	37	Fr2	V	V	V	V	V	V
	38	38	38	Fr2	R	R	R	R	R	R
	39	39	39	Fr2	Q	Q	Q	Q	Q	Q
50	40	40	40	Fr2	T	T	T	T	T	T
	41	41	41	Fr2	P	P	P	P	P	P
	42	42	42	Fr2	E	E	D	E	E	E

55

**Table 13**  
**Humanized 14G8 Vh Regions**

	Chelonia residue #	Kabat residue #	Linear residue #	FR or CDR	Mouse 14G8 VH (SEQ ID NO: 61)	Hu VH Acceptor Fr Ace# AAD30410.1 (SEQ ID NO: 63)	Hu VH Acceptor Fr Ace# AAK82494.1 (SEQ ID NO: 4)	Hu14G8 VHv1 (SEQ ID NO: 64)	Hu14G8 VHv2 (SEQ ID NO: 65)	Hu14G8 VHv3 (SEQ ID NO: 66)
5										
10										
15										
20										
25										
30										
35										
40										
45										
50										
55										
	43	43	43	Fr2	K	K	K	K	K	K
	44	44	44	Fr2	R	R	R	R	R	R
	45	45	45	Fr2	L	L	L	L	L	L
	46	46	46	Fr2	E	E	E	E	E	E
	47	47	47	Fr2	L	W	W	L	L	L
	48	48	48	Fr2	V	V	V	V	V	V
	49	49	49	Fr2	A	A	A	A	A	A
	50	50	50	CDR-H2	E	A	T	E	E	E
	51	51	51	CDR-H2	I	I	I	I	I	I
	52	52	52	CDR-H2	N	S	S	N	N	N
	52A	52A	53	CDR-H2	N	T	S	N	N	N
	53	53	54	CDR-H2	S	D	G	S	S	S
	54	54	55	CDR-H2	G	G	G	G	G	G
	55	55	56	CDR-H2	D	S	S	D	D	D
	56	56	57	CDR-H2	T	F	Y	T	T	T
	57	57	58	CDR-H2	T	I	I	T	T	T
	58	58	59	CDR-H2	Y	Y	Y	Y	Y	Y
	59	59	60	CDR-H2	Y	Y	Y	Y	Y	Y
	60	60	61	CDR-H2	P	A	P	P	P	P
	61	61	62	CDR-H2	D	D	D	D	D	D
	62	62	63	CDR-H2	T	T	S	T	T	T
	63	63	64	CDR-H2	V	V	V	V	V	V
	64	64	65	CDR-H2	K	K	K	K	K	K
	65	65	66	CDR-H2	G	G	G	G	G	G
	66	66	67	Fr3	R	R	R	R	R	R
	67	67	68	Fr3	F	F	F	F	F	F
	68	68	69	Fr3	T	T	T	T	T	T
	69	69	70	Fr3	I	I	I	I	I	I
	70	70	71	Fr3	S	S	S	S	S	S
	71	71	72	Fr3	R	G	R	R	R	R

55

**Table 13**  
**Humanized 14G8 Vh Regions**

Chelvia residue #	Kabat residue #	Linear residue #	FR or CDR	Murine 14G8 VH (SEQ ID NO: 61)	Hu VH Acceptor Fr Ace# AAD30410.1 (SEQ ID NO: 63)	Hu VH Acceptor Fr Ace# AAK82494.1 (SEQ ID NO: 4)	Hu14G8 VHV1 (SEQ ID NO: 64)	Hu14G8 VHV2 (SEQ ID NO: 65)	Hu14G8 VHV3 (SEQ ID NO: 66)
72	72	73	Fr3	D	D	D	D	D	D
73	73	74	Fr3	N	N	N	N	N	N
74	74	75	Fr3	A	S	A	A	A	A
75	75	76	Fr3	K	K	K	K	K	K
76	76	77	Fr3	N	N	N	N	N	N
77	77	78	Fr3	T	T	T	T	T	T
78	78	79	Fr3	L	L	L	L	L	L
79	79	80	Fr3	Y	Y	Y	Y	Y	Y
80	80	81	Fr3	L	L	L	L	L	L
81	81	82	Fr3	Q	Q	Q	Q	Q	Q
82	82	83	Fr3	M	M	M	M	M	M
82A	82A	84	Fr3	S	N	S	S	S	N
82B	82B	85	Fr3	S	S	S	S	S	S
82C	82C	86	Fr3	L	L	L	L	L	L
83	83	87	Fr3	K	R	K	K	K	R
84	84	88	Fr3	S	A	S	S	S	A
85	85	89	Fr3	E	E	E	E	E	E
86	86	90	Fr3	D	D	D	D	D	D
87	87	91	Fr3	T	T	T	T	T	T
88	88	92	Fr3	A	A	A	A	A	A
89	89	93	Fr3	M	V	M	M	M	V
90	90	94	Fr3	Y	Y	Y	Y	Y	Y
91	91	95	Fr3	Y	Y	Y	Y	Y	Y
92	92	96	Fr3	C	C	C	C	C	C
93	93	97	Fr3	A	A	A	A	A	A
94	94	98	Fr3	R	K	R	R	R	R
95	95	99	CDR-H3	H	D	L	H	H	H
96	96	100	CDR-H3	Y	R	Y	Y	Y	Y
97	97	101	CDR-H3	Y	G	Y	Y	Y	Y
98	98	102	CDR-H3	Y	I	G	Y	Y	Y

**Table 13**  
**Humanized 14G8 Vh Regions**

Chordata residue #	Kabat residue #	Linear residue #	FR of CDR	Mouse 14G8 VH (SEQ ID NO: 61)	Hu VH Acceptor Fr Acc# AAD30410.1 (SEQ ID NO: 63)	Hu VH Acceptor Fr Acc# AAX82494.1 (SEQ ID NO: 4)	Hu14G8 VHV1 (SEQ ID NO: 64)	Hu14G8 VHV2 (SEQ ID NO: 65)	Hu14G8 VHV3 (SEQ ID NO: 66)
99	99	103	CDR-H3	G	D	Y	G	G	G
100	100	104	CDR-H3	G	A	R	G	G	G
100A	100A	105	CDR-H3	G	T	Y	G	G	G
100B	100B	106	CDR-H3	Y	A	Y	Y	Y	Y
100C	100C	107	CDR-H3	G	Q	F	G	G	G
100D	100D	108	CDR-H3	G	V	-	G	G	G
100E	100E	109	CDR-H3	W	G	-	W	W	W
100F	100F	110	CDR-H3	F	R	-	F	F	F
100G	100G	111	CDR-H3	F	F	-	F	F	F
101	101	112	CDR-H3	D	D	D	D	D	D
102	102	113	CDR-H3	V	P	Y	V	V	V
103	103	114	Fr4	W	W	W	W	W	W
104	104	115	Fr4	G	G	G	G	G	G
105	105	116	Fr4	T	Q	Q	T	Q	Q
106	106	117	Fr4	G	G	G	G	G	G
107	107	118	Fr4	T	T	T	T	T	T
108	108	119	Fr4	T	L	M	L	L	L
109	109	120	Fr4	V	V	V	V	V	V
110	110	121	Fr4	T	T	T	T	T	T
111	111	122	Fr4	V	V	V	V	V	V
112	112	123	Fr4	S	S	S	S	S	S
113	113	124	Fr4	S	S	S	S	S	S

**Table 14**  
**Humanized 14G8 Vh Regions**

Chenbala residue #	Kabat residue #	Linear residue #	FR of CDR	Mouse 14G8 VL (SEQ ID NO: 70)	Hu VL Acceptor Fr Acc# ABA71374.1 (SEQ ID NO: 72)	Hu VL Acceptor Fr Acc# ABC66952.1 (SEQ ID NO: 73)	Hu14G8 VLv1 (SEQ ID NO: 74)	Hu14G8 VLv2 (SEQ ID NO: 75)	Hu14G8 VLv3 (SEQ ID NO: 76)
1	1	1	FR1	D	D	D	D	D	D
2	2	2	FR1	F	F	F	F	F	F
3	3	3	FR1	V	V	V	V	V	V
4	4	4	FR1	M	M	M	M	M	M
5	5	5	FR1	T	T	T	T	T	T
6	6	6	FR1	Q	Q	Q	Q	Q	Q
7	7	7	FR1	A	T	S	S	S	S
8	8	8	FR1	A	P	P	A	P	P
9	9	9	FR1	P	L	L	P	L	L
10	10	10	FR1	S	S	S	S	S	S
11	11	11	FR1	V	L	L	L	L	L
12	12	12	FR1	P	P	P	P	P	P
13	13	13	FR1	V	V	V	V	V	V
14	14	14	FR1	T	T	T	T	T	T
15	15	15	FR1	P	P	P	P	P	P
16	16	16	FR1	G	G	G	G	G	G
17	17	17	FR1	E	E	E	E	E	E
18	18	18	FR1	S	S	P	S	P	P
19	19	19	FR1	V	A	A	V	A	A
20	20	20	FR1	S	S	S	S	S	S
21	21	21	FR1	F	F	F	F	F	F
22	22	22	FR1	S	S	S	S	S	S
23	23	23	FR1	C	C	C	C	C	C
24	24	24	CDR-L1	R	R	R	R	R	R
25	25	25	CDR-L1	S	S	S	S	S	S
26	26	26	CDR-L1	N	S	N	N	N	S
27	27	27	CDR-L1	K	Q	Q	K	P	P
27A	27A	28	CDR-L1	S	S	S	S	S	S

55

**Table 14**  
**Humanized 14G8 Vk Regions**

Chothia residue #	Kabat residue #	Linear residue #	FR or CDR	Mouse 14G8 VL (SEQ ID NO: 70)	Hu VL Acceptor Fr Acc# ABA71374.1 (SEQ ID NO: 72)	Hu VL Acceptor Fr Acc# ABC68952.1 (SEQ ID NO: 73)	Hu14G8 VLv1 (SEQ ID NO: 74)	Hu14G8 VLv2 (SEQ ID NO: 75)	Hu14G8 VLv3 (SEQ ID NO: 76)
5									
10									
15									
20									
25									
30									
35									
40									
45									
50									
55									
	27B	29	CDR-L1	L	L	L	L	L	L
	27C	30	CDR-L1	L	R	L	L	L	L
	27D	31	CDR-L1	H	H	Y	H	H	H
	27E	32	CDR-L1	S	Y	S	S	S	S
	27F	33	CDR-L1	N	S	N	N	N	N
	28	34	CDR-L1	G	G	G	G	G	G
	29	35	CDR-L1	N	Y	Y	N	N	N
	30	36	CDR-L1	I	I	N	I	I	I
	31	37	CDR-L1	Y	Y	Y	Y	Y	Y
	32	38	CDR-L1	L	I	L	L	L	L
	33	39	CDR-L1	Y	D	D	Y	Y	Y
	34		CDR-L1						
	35	40	Fr2	W	W	W	W	W	W
	36	41	Fr2	F	Y	Y	F	F	F
	37	42	Fr2	L	L	L	L	L	L
	38	43	Fr2	Q	Q	Q	Q	Q	Q
	39	44	Fr2	R	K	K	K	K	K
	40	45	Fr2	P	P	P	P	P	P
	41	46	Fr2	G	G	G	G	G	G
	42	47	Fr2	Q	Q	Q	Q	Q	Q
	43	48	Fr2	S	S	S	S	S	S
	44	49	Fr2	P	P	P	P	P	P
	45	50	Fr2	Q	Q	Q	Q	Q	Q
	46	51	Fr2	L	V	L	L	L	L
	47	52	Fr2	L	L	L	L	L	L
	48	53	Fr2	I	I	I	I	I	I
	49	54	Fr2	Y	Y	Y	Y	Y	Y
	50	55	CDR-L2	R	L	S	R	R	R
	51	56	CDR-L2	V	G	G	V	V	V
	52	57	CDR-L2	S	S	S	S	S	S

**Table 14**  
**Humanized 14G8 Vk Regions**

5

10

15

20

25

30

35

40

45

50

55

Chothia residue #	Kabat residue #	Linear residue #	FR or CDR	Mouse 14G8 VL (SEQ ID NO: 70)	Hu VL Acceptor Fr Acc# ABA71374.1 (SEQ ID NO: 72)	Hu VL Acceptor Fr Acc# ABC68952.1 (SEQ ID NO: 73)	Hu14G8 VLv1 (SEQ ID NO: 74)	Hu14G8 VLv2 (SEQ ID NO: 75)	Hu14G8 VLv3 (SEQ ID NO: 76)
53	53	58	CDR-L2	N	N	N	N	N	N
54	54	59	CDR-L2	L	R	R	L	L	L
55	55	60	CDR-L2	A	A	A	A	A	A
56	56	61	CDR-L2	S	S	S	S	S	S
57	57	62	Fr3	G	G	G	G	G	G
58	58	63	Fr3	V	V	V	V	V	V
59	59	64	Fr3	P	P	P	P	P	P
60	60	65	Fr3	D	D	D	D	D	S
61	61	66	Fr3	R	R	R	R	R	R
62	62	67	Fr3	F	F	F	F	F	F
63	63	68	Fr3	S	S	S	S	S	S
64	64	69	Fr3	G	G	G	G	G	G
65	65	70	Fr3	S	S	S	S	S	S
66	66	71	Fr3	G	G	G	G	G	G
67	67	72	Fr3	S	S	S	S	S	S
68	68	73	Fr3	G	G	G	G	G	G
69	69	74	Fr3	T	T	T	T	T	T
70	70	75	Fr3	A	D	D	A	D	D
71	71	76	Fr3	F	F	F	F	F	F
72	72	77	Fr3	T	T	T	T	T	T
73	73	78	Fr3	L	L	L	L	L	L
74	74	79	Fr3	R	K	K	K	K	K
75	75	80	Fr3	I	I	I	I	I	I
76	76	81	Fr3	S	S	S	S	S	S
77	77	82	Fr3	R	R	R	R	R	R
78	78	83	Fr3	V	V	V	V	V	V
79	79	84	Fr3	E	E	E	E	E	E
80	80	85	Fr3	A	A	A	A	A	A
81	81	86	Fr3	E	E	E	E	E	E
82	82	87	Fr3	D	D	D	D	D	D

**Table 14**  
**Humanized 14G8 VL Regions**

5

10

15

20

25

30

35

40

45

50

55

Chothia residue #	Kabat residue #	Linear residue #	FR or CDR	Mouse 14G8 VL (SEQ ID NO: 70)	Hu VL Acceptor Fr Acc# ABA71374.1 (SEQ ID NO: 72)	Hu VL Acceptor Fr Acc# ABC66952.1 (SEQ ID NO: 73)	Hu14G8 VLv1 (SEQ ID NO: 74)	Hu14G8 VLv2 (SEQ ID NO: 75)	Hu14G8 VLv3 (SEQ ID NO: 76)
83	83	88	Fc3	V	V	V	V	V	V
84	84	89	Fc3	G	G	G	G	G	G
85	85	90	Fc3	V	V	V	V	V	V
86	86	91	Fc3	Y	Y	Y	Y	Y	Y
87	87	92	Fc3	Y	Y	Y	Y	Y	Y
88	88	93	Fc3	C	C	C	C	C	C
89	89	94	CDR-L3	M	M	M	M	M	M
90	90	95	CDR-L3	Q	Q	Q	Q	Q	Q
91	91	96	CDR-L3	H	A	A	H	H	H
92	92	97	CDR-L3	L	L	L	L	L	L
93	93	98	CDR-L3	E	Q	Q	E	E	E
94	94	99	CDR-L3	Y	T	S	Y	Y	Y
95	95	100	CDR-L3	P	V	P	P	P	P
96	96	101	CDR-L3	L	C	Y	L	L	L
97	97	102	CDR-L3	T	S	T	T	T	T
98	98	103	Fc4	F	F	F	F	F	F
99	99	104	Fc4	G	G	G	G	G	G
100	100	105	Fc4	A	Q	Q	Q	Q	Q
101	101	106	Fc4	G	G	G	G	G	G
102	102	107	Fc4	T	T	T	T	T	T
103	103	108	Fc4	K	K	K	K	K	K
104	104	109	Fc4	L	L	L	L	L	L
105	105	110	Fc4	E	E	E	E	E	E
106	106	111	Fc4	L	I	I	I	I	I
106A	106A	112	Fc4	K	K	K	K	K	K
107	107	113	Fc4	R	R	R	R	R	R

Table 15

<b>V<sub>H</sub>, V<sub>L</sub> Backmutations and Other Mutations</b>			
<b>V<sub>H</sub> or V<sub>L</sub> Variant</b>	<b>V<sub>H</sub> or V<sub>L</sub> Exon Acceptor Sequence</b>	<b>Donor Framework Residues</b>	<b>Kabat CDR Residues</b>
Hu14G8VHv1 (SEQ ID NO: 64)	NCBI accession codes AAD30410.1 and AAX82494.1 (SEQ ID NOS: 63 and 4, respectively)	H1, H3, H47, H105	-
Hu14G8VHv2 (SEQ ID NO: 65)	NCBI accession codes AAD30410.1 and AAX82494.1 (SEQ ID NOS: 63 and 4, respectively)	H1, H47	-
Hu14G8VHv3 (SEQ ID NO: 66)	NCBI accession codes AAD30410.1 and AAX82494.1 (SEQ ID NOS: 63 and 4, respectively)	H1, H47	-
Hu14G8VLv1 (SEQ ID NO: 74)	NCBI accession codes ABA71374.1 and ABC66952.1 (SEQ ID NOS: 72 and 73, respectively)	L8, L9, L19, L36, L70	-
Hu14G8VLv2 (SEQ ID NO: 75)	NCBI accession codes ABA71374.1 and ABC66952.1 (SEQ ID NOS: 72 and 73, respectively)	L36	-
Hu14G8VLv3 (SEQ ID NO: 76)	NCBI accession codes ABA71374.1 and ABC66952.1 (SEQ ID NOS: 72 and 73, respectively)	L36, L60	L26

Table 16

<b>Kabat Numbering of Framework and Kabat CDR Residues for Backmutations and Other Mutations in Humanized 14G8 Antibody VH Regions</b>						
<b>Residue</b>	<b>AAD30410.1 Heavy Chain</b>	<b>AAX82494.1 Heavy Chain</b>	<b>Mouse 14G8</b>	<b>Hu14G8 VHv1</b>	<b>Hu14G8 VHv2</b>	<b>Hu14G8 VHv3</b>
H1	Q	Q	E	E	E	E
H3	Q	Q	K	K	Q	Q
H47	W	W	L	L	L	L
H105	Q	Q	T	T	Q	Q

Table 17

<b>Kabat Numbering of Framework and Kabat CDR Residues for Backmutations and Other Mutations in Humanized 14G8 Antibody VL Regions</b>						
<b>Residue</b>	<b>ABA71374.1 Light Chain</b>	<b>ABC66952.1 Light Chain</b>	<b>Mouse 14G8</b>	<b>Hu14G8 VLv1</b>	<b>Hu14G8 VLv2</b>	<b>Hu14G8 VLv3</b>
L8	P	P	A	A	P	P
L9	L	L	P	P	L	L
L19	A	A	V	V	A	A
L26	S	N	N	N	N	S
L36	Y	Y	F	F	F	F
L60	D	D	D	D	D	s
L70	D	D	A	A	D	D

**[0280]** An alignment of the murine 14G8 Vh sequence (SEQ ID NO: 61) with the mouse model sequence (1MQK H; SEQ ID NO: 62), the human acceptor sequences (AAD30410.1 and AAX82494.1; SEQ ID NOS: 63 and 4, respectively), and the Hu14G8VHv1, Hu14G8VHv2, and Hu14G8VHv3 sequences (SEQ ID NOS: 64-66, respectively), is shown in Figure 2A. The CDRs as defined by Kabat are enclosed in boxes, except that the first enclosed box is a composite of the Chothia CDR-H1 and the Kabat CDR-H1, with the Kabat CDR-H1 underlined and bolded. Positions at which canonical, vernier, or interface residues differ between mouse and human acceptor sequences are candidates for substitution. Examples of vernier/CDR foundation residues include residues 2, 49, 69, 71, 75, 80, and 94 by Kabat numbering in Table 13. Examples of canonical/CDR interacting residues include residues 24, 48, and 73 by Kabat numbering in Table 13. Examples of interface/packing (VH+VL) residues include residues 37, 39, 44, 47, 91, 93, and 103 by Kabat numbering in Table 13.

**[0281]** An alignment of the murine 9D5 Vk 14G8 (SEQ ID NO: 70) with the mouse model sequence (1MJU L; SEQ ID NO: 71), the human acceptor sequences (ABA71374.1 and ABC66952.1; SEQ ID NOS: 72 and 73, respectively), and the Hu14G8VLv1, Hu14G8VLv2, Hu14G8VLv3 sequences (SEQ ID NOS: 74-76, respectively), is shown in Figure 2B. The CDRs as defined by Kabat are enclosed in boxes, except that the first enclosed box is a composite of the Chothia CDR-H1 and the Kabat CDR-H1, with the Kabat CDR-H1 underlined and bolded. Positions at which canonical, vernier, or interface residues differ between mouse and human acceptor sequences are candidates for substitution. Examples of vernier/CDR foundation residues include residues 4, 35, 46, 49, 66, 68, and 69 by Kabat numbering in Table 14. Examples of canonical/CDR interacting residues include residues 2, 48, 64, and 71 by Kabat numbering in Table 14. Examples of interface/packing (VH+VL) residues include residues 36, 38, 44, 47, 87, and 98 by Kabat numbering in Table 14.

**[0282]** The rationales for selection of the positions indicated in Tables 15 and 17 in the light chain variable region as candidates for substitution are as follows.

**[0283]** P8A: Proline is critical to structure conformation, and loss or gain of proline may affect the structure. Backmutations were designed to avoid "gain of proline." However, since A at this position is rare in human IgG frameworks, P was included in some versions.

**[0284]** L9P: Proline is critical to structure conformation, and loss or gain of proline may affect the structure. Backmutations were designed to avoid "loss of proline." However, since P at this position is rare in human IgG frameworks, L was included in some versions.

**[0285]** A19V: A and V have similar frequencies in human frameworks. Both residues were included in separate versions.

**[0286]** N26S: L26 represents an N-glycosylation site in light chain CDR1. Mutation to S reduces N-glycosylation and yields a more heterogeneous product.

**[0287]** Y36F: This is an interface residue. Substitution was made in some versions to keep the murine residue at this interface residue.

**[0288]** D60S: D and S are similar in frequency in human IgG frameworks. S is applied by most commercially available therapeutic antibodies, so D was replaced with S in some versions.

**[0289]** D70A: D is much more frequent than A in human IgG frameworks. However, D will form ionic bond with R24 in light chain CDR1, so both A and D were included in separate versions.

**[0290]** The rationales for selection of the positions indicated in Tables 15 and 16 in the heavy chain variable region as candidates for substitution are as follows.

**[0291]** Q1E: Glutamate (E) conversion to pyroglutamate (pE) occurs more slowly than from glutamine (Q). Because of the loss of a primary amine in the glutamine to pE conversion, antibodies become more acidic. Incomplete conversion produces heterogeneity in the antibody that can be observed as multiple peaks using charge-based analytical methods. Heterogeneity differences may indicate a lack of process control.

**[0292]** Q3K: Q is more frequent in human IgG frameworks. K is less frequent but not rare. Both residues were included in separate versions.

**[0293]** W47L: This is an interface residue. Substitution was made in some versions to keep the murine residue at this interface residue.

**[0294]** Q105T: T may form a hydrogen bond with K3, so T was included in some versions to maintain the conformational structure of VH.

**[0295]** Because two human acceptor frameworks (ABA71374.1 and ABC66952.1 (SEQ ID NOS: 72 and 73, respectively)) were used for humanization of the 14G8 light chain variable region, there are certain framework positions that have different amino acids in the two acceptor frameworks. Humanized versions of the 14G8 light chain variable region can include either amino acid at these residues. An example of a residue in which these two acceptors differ is position L18 (S or P) by Kabat numbering. The rationale for choosing one amino acid or the other at this position is as follows.

**[0296]** L18 (S or P): P at this position is less frequent. Mutation to S was tried to alleviate loop distortion.

**[0297]** Because two human acceptor frameworks (AAD30410.1 and AAX82494.1 (SEQ ID NOS: 63 and 4, respectively)) were also used as acceptor sequences for humanization of the 14G8 mature heavy chain variable region, there are certain framework positions that have different amino acids in the two acceptor frameworks. Humanized versions

of the 14G8 heavy chain variable region can include either amino acid at these residues. Examples of residues in which these two acceptors differ include positions H82a (N or S), H83 (R or K), H84 (A or S), and H89 (V or M) by Kabat numbering. The rationales for choosing one amino acid or the other at these positions are as follows.

**[0298]** H82a (N or S): N and S have similar frequency on 82a. They are also human residues in two human VH templates. S and N were included in separate versions.

**[0299]** H83 (R or K): K in framework region 3 is close to the CDR binding surface area. Replacing it with R might obstruct placement of antigen. R is most frequent and K is third frequent in human frameworks at this position. R and K were included in separate versions.

**[0300]** H84 (A or S): S and A are present in the two frameworks considered, so each residue was tried in separate versions.

**[0301]** H89 (V or M): V and M are present in the two frameworks considered, so each was tried in separate versions.

**[0302]** The three humanized light chain variable region variants and three humanized heavy chain variable region variants are as follows:

**[0303]** Hu14G8VL version 1 (P8A, L9P, A19V, Y36F, D70A substitutions in lowercase):

DIVMTQ<sup>Sap</sup>SLPVT<sup>PGES</sup>vSISCRSNKSL<sup>LH</sup>SN<sup>NTY</sup>LYW<sup>LQ</sup>KPGQSPQL<sup>LI</sup>YRVSNLASG  
VPDRFSGSGSGT<sup>FT</sup>LKISRVEAEDVGVYYCMQHLEYP<sup>LT</sup>FGQG<sup>TK</sup>LEIKR (SEQ ID NO:  
74)

**[0304]** Hu14G8VL version 2 (L36F substitution in lowercase):

DIVMTQSP<sup>LS</sup>LPVT<sup>PGEP</sup>ASISCRSNKSL<sup>LH</sup>SN<sup>NTY</sup>LYW<sup>LQ</sup>KPGQSPQL<sup>LI</sup>YRVSNLAS  
GVPDRFSGSGSGT<sup>DF</sup>TLKISRVEAEDVGVYYCMQHLEYP<sup>LT</sup>FGQG<sup>TK</sup>LEIKR (SEQ ID  
NO: 75)

**[0305]** Hu14G8VL version 3 (N26S, L36F, and D60S substitutions in lowercase):

DIVMTQSP<sup>LS</sup>LPVT<sup>PGEP</sup>ASISCRS<sup>S</sup>KSL<sup>LH</sup>SN<sup>NTY</sup>LYW<sup>LQ</sup>KPGQSPQL<sup>LI</sup>YRVSNLAS  
GVP<sup>S</sup>RFSGSGSGT<sup>DF</sup>TLKISRVEAEDVGVYYCMQHLEYP<sup>LT</sup>FGQG<sup>TK</sup>LEIKR (SEQ ID  
NO: 76)

**[0306]** Hu14G8VH version 1 (Q1E, Q3K, W47L, and Q105T substitutions in lowercase):

eV<sup>L</sup>LVESGGGLVQPGGSLK<sup>L</sup>SCAASGFTFSSYTMSWVRQTPEKRLEIVAEINNSGDTTYY  
PDTVKGRFTISRDN<sup>AKNT</sup>LYLQ<sup>M</sup>SSLK<sup>SED</sup>TAMYYCARHYYYGGGYGGWFFD<sup>V</sup>WG<sup>I</sup>G  
TLVTVSS (SEQ ID NO: 64)

**[0307]** Hu14G8VH version 2 (Q1E and W47L substitutions in lowercase):

eV<sup>Q</sup>LVESGGGLVQPGGSLK<sup>L</sup>SCAASGFTFSSYTMSWVRQTPEKRLEIVAEINNSGDTTYY  
PDTVKGRFTISRDN<sup>AKNT</sup>LYLQ<sup>M</sup>SSLK<sup>SED</sup>TAMYYCARHYYYGGGYGGWFFD<sup>V</sup>WG<sup>Q</sup>  
GTLVTVSS (SEQ ID NO: 65)

**[0308]** Hu14G8VH version 3 (Q1E and W47L substitutions in lowercase):

eV<sup>Q</sup>LVESGGGLVQPGGSLK<sup>L</sup>SCAASGFTFSSYTMSWVRQTPEKRLEIVAEINNSGDTTYY  
PDTVKGRFTISRDN<sup>AKNT</sup>LYLQ<sup>M</sup>NSLRAEDTAVYYCARHYYYGGGYGGWFFD<sup>V</sup>WG<sup>Q</sup>  
GTLVTVSS (SEQ ID NO: 66)

#### Example 10. Binding Kinetic Analysis of Humanized 14G8 Antibodies

**[0309]** Binding kinetics of three humanized 14G8 variants and murine 14G8 to recombinant human TTR F87M/L110M

were characterized by Biacore, as shown in Table 18. Anti-human (GE Healthcare) was immobilized on sensor chip C5 (lacking dextran chains) via amine coupling following the instructions provided in the GE Healthcare anti-human kit, and mAbs were captured to a level to ensure a maximum binding of analyte of 30-50 RU. Various concentrations of analyte (recombinant human TTR F87M/L110M) were passed over the captured ligand at 50  $\mu$ l/min in running buffer (HBS + 0.05% P-20, 1 mg/mL BSA) in 3-fold dilutions. For each concentration, the reaction proceeded for a time frame allowing for the higher analyte concentrations to reach equilibrium during association, as well as at least 10% of signal to decay during dissociation. At least one concentration (not the highest or lowest) was run in duplicate. Concentration ranges of analyte were selected based on preliminary experimentation to span at least 10-fold above  $K_D$  to 10-fold below  $K_D$ .

**[0310]** Table 18 summarizes the association rate ( $k_a$ ), dissociation rate ( $k_d$ ), and binding affinity constant ( $K_D$ ) of the mouse 14G8, chimeric-14G8, Hu14G8 H2L1, Hu14G8 H2L2, Hu14G8H2L3, and Hu14G8 H3L1 for recombinant human TTR F87M/L110M. As shown in Table 18, Hu14G8 antibodies and mouse 14G8 have similar TTR binding affinities.

Table 18

Antigen Binding Affinity of 14G8 Antibodies to Human TTR (F87M/L110M)				
Antibody	$k_a$ (1/Ms)	$k_d$ (1/s)	$K_D$ (M)	Rmax
Murine 14G8	2.88E+04	5.36E-04	1.86E-08	27.13
Chimeric 14G8	2.51E+05	3.46E-04	1.38E-09	51.06
Hu14G8 H2L1	3.07E+05	6.93E-04	2.26E-09	42.03
Hu14G8 H2L2	3.26E+05	6.83E-04	2.10E-09	51.43
Hu14G8 H2L3	2.52E+05	3.18E-04	1.24E-09	42.41
Hu14G8 H3L1	3.86E+05	6.77E-04	1.76E-09	44.05

**[0311]** Murine 14G8, chimeric 14G8, and other humanized version of 14G8 with no changes in LCDR1 showed additional light chains when run on SDS-PAGE under reducing conditions. Sequence analysis revealed that there is an N-glycosylation site in LCDR1 at residue 26 by Kabat numbering. This N-glycosylation site may cause potential heterogeneity problems during manufacture. Mutation of the N at residue L26 in Hu14G8VLv2 yielded Hu14G8VLv3. Humanized antibodies having Hu14G8VLv3 show only one light chain when run on SDS-PAGE under reducing conditions, thereby eliminating the potential heterogeneity problem. The antigen binding affinity of Hu14G8 H2L3 is similar to that of the murine parent antibody and the chimeric antibody.

### Example 11. Materials and Methods

#### a. Antibody Generation Protocol

**[0312]** Mice were immunized weekly with the antigenic peptides TTR-MAP, TTR89-97-N-KLH or TTR89-97-C-KLH in RIBI adjuvant or monthly in TiterMax adjuvant. Three to four days prior to fusion, selected mice were given a final IV boost with immunogen in saline solution. Spleen were homogenized to prepare splenocytes and fused with SP2/0 myeloma cells using a standard electrofusion protocol. Fused cells in selection media were plated in 96-well plates and screened after 7-10 days.

#### b. Antibody Screening Protocol

**[0313]** Hybridoma selection was based on the following ELISA screen: 96-well ELISA plates were coated with chicken anti-His, 1  $\mu$ g/mL PBS and incubated for 1 hour. Plates were blocked with of 1% BSA/PBS solution, 200  $\mu$ L/well for 15 minutes then 0.5  $\mu$ g/mL pH4-TTR, 50  $\mu$ L/well was added and incubated for 1 hour. pH4-TTR is TTR that has been subjected to low pH (50 mM sodium acetate, pH 4.0) in order to dissociate/aggregate TTR, exposing the TTR89-97 epitope. Plates were washed twice with TBS-T. Supernatant from fusion plates was added, 50  $\mu$ L/well and incubated for 1 hour. Plates were washed twice with TBS-T. The detection antibody, goat anti-mouse (IgG1, 2a, 2b, 3 specific)-HRP diluted 1:5,000 in 0.5% BSA/PBS/TBS-T, 50  $\mu$ L/well was added and incubated for 1 hour. Finally, plates were washed five times with TBS-T and TMB substrate, 100  $\mu$ L/well was added. After 15 minutes, substrate development was stopped with 2N Sulfuric Acid, 50  $\mu$ L/well. Plates were read at 450 nm. Wells with an O.D. > 1.0 were selected and cells were transferred to a 24-well plate. After 3 days of growth, clones were counter screened with the above assay to confirm binding, and substituting native TTR for pH4-TTR as a negative counter screen, allowing for selection of clones producing TTR mAbs specific for non-native forms of TTR.

### c. Antibody Expression Protocols

**[0314]** CMV driven light chain and heavy chain plasmids carrying humanized monoclonal antibody sequences were transfected into CHO-S1 cells (Life Technology). Dual selection was applied to make a selected pool. Conditioned media was assayed for titer, binding and analyzed by SDS-PAGE/Western blotting. Selected pools were used for clone generation using Clonepix system (Molecular Devices). Clones were ranked based on antibody titer. Selected clones were expanded and banked.

**[0315]** The highest producing clone was expanded in shake flasks and the culture was used to inoculate 10-25L Wave bag cultures. A mixture of FreeStyle-CHO, CD OptiCHO and FreeStyle F17 expression media supplemented with Glutamax (media and Glutamax from Life Technology) was used for shake flask as well as for Wave bag cultures. Batch culture was made using a Wave Bioreactor (GE Healthcare) at 37°C, 7% CO<sub>2</sub> under constant agitation. Samples were drawn periodically to monitor cell number, viability and antibody production. Supplementation with Cell Boost (HyClone) was made if needed. The batch culture was harvested when cell viability starts to decline below 90% (5 -7 days).

### d. Antibody Purification Protocol

**[0316]** The cell culture was harvested after first allowing the cells in suspension to settle down to the bottom of the Wave bag via gravity at 4°C. Harvested media was clarified through a depth filter (Millistak Pod COHC, Millipore), concentrated 10-fold by tangential flow filtration (Pelicon 2PLC 30K, Millipore) and sterile filtered through a 0.2 µm filter (Opticap XL, Millipore). The concentrated conditioned media was then loaded onto a Protein G Sepharose Fast Flow column (GE Lifesciences) pre-equilibrated in 1xPBS, pH 7.4 using an FPLC (Akta Avant, GE Lifesciences). Unbound proteins were washed off the column with 5-10 column volumes of 1xPBS, pH 7.4 until the OD<sub>280</sub> reached baseline. The bound antibody was eluted from the column with 2 column volumes of IgG Elution Buffer (Thermo Scientific). Elution fractions were collected and pH neutralized with 2M Tris, pH 9.0 (60 µL per 1ml elution).

**[0317]** Antibody-containing fractions were pooled and dialyzed overnight at 4°C against 1xPBS, pH 7.4. The dialyzed sample was then sterilized by ultrafiltration through a 0.2 µm PES filter and stored at 4°C. The final protein concentration was determined by bicinchoninic acid (BCA) using bovine gamma-globulin as the protein standard (Thermo Scientific).

### e. Recombinant TTR Expression and Purification Protocols

**[0318]** E. coli (BL21-A1) cells were transformed with a pET21a(+) plasmid containing a TTR insert (Met-hTTR-(His)<sub>6</sub> or a TTR variant containing an F87M/L110M double mutation. Cells were grown in 2YT broth containing 100 µg/ml ampicillin. Expression of TTR was induced overnight at 20°C in the presence of 1mM IPTG and 0.05% arabinose.

**[0319]** The cells were collected by centrifugation at 4000 x g for 10 min. and stored at -80°C until used. 10-15g cell pellets were thawed and lysed in 50ml Buffer A (1xPBS containing 500 mM NaCl, 20mM imidazole) by processing through an LV-1 high-shear processor (Microfluidics, Inc.). Lysed cells were centrifuged at 12,000 x g for 15 min, filtered through a 0.2 µm PES filter prior to purification on a His-Trap HP column (GE Lifesciences). After loading, the column was washed with 10 c.v. of Buffer A and eluted with Buffer B (1xPBS with 500 mM NaCl, 500 mM imidazole). Peak fractions corresponding to TTR were collected, dialyzed against 1xPBS and stored at -80°C until used.

### f. TTR Antigen Preparation

**[0320]** Native TTR antigen was prepared by diluting a concentrated stock of recombinant TTR-6His to a final concentration of 2.5 µg/ml in 1x PBS. pH4-treated TTR was generated by incubating recombinant TTR at a concentration of 0.2 mg/ml in 50mM sodium acetate, pH 3.95 for 72 hours at room temperature. Under these conditions, TTR dissociates into mixture of TTR monomers and aggregated forms that are structurally distinct from native TTR. The pH4-TTR was then diluted to a final concentration of 2.5 µg/ml in 1xPBS immediately before use in the assay. 96-well plates (Costar #3690) were coated at room temperature with 50 µl per well of 1.0 µg/ml chicken-anti-his polyclonal antibody (Abcam #Ab9107) in 1xPBS for 1 hr. The coating solution was discarded and the plate was blocked with a 250 µl per well volume of 1x BSA-containing block buffer diluted in 1xPBS (G-Biosciences #786-193) for 1hr.

### g. ELISA Protocol

**[0321]** Coated and blocked 96-well plates were treated with 50 µl per well of 2.5 µg/ml TTR antigen (either native TTR or pH4-TTR) for 1hr. at room temperature. The plates were then washed two times with 250 µl per well of wash buffer (1x Tris Buffered Saline containing 0.05% Tween-20). Washed plates were then treated with 50 µl per well of the appropriate anti-TTR monoclonal antibody at concentrations ranging from 0.31 to 2.5 µg/ml, for 1hr.

**[0322]** The treated plates were washed 3 times with 250 µl per well wash buffer. After washing, the plates were treated

for 1 hr. with 50  $\mu$ l per well of detection antibody comprising a 1:5,000 dilution of peroxide-conjugated goat-anti-mouse (Jackson ImmunoResearch #115-035-164) in 1xPBS. The plate was then washed 3 times prior to the addition of 100  $\mu$ l per well TMB substrate (Rockland). The HRP reaction was allowed to proceed at room temperature for 15min. before quenching with a 50  $\mu$ l per well volume of IN H<sub>2</sub>SO<sub>4</sub>. Spectroscopic absorbance was measured at a wavelength of 450nm.

5

#### *h. SDS-PAGE*

**[0323]** Electrophoresis on SDS-polyacrylamide gels was carried out as follows. 0.1-1  $\mu$ g TTR or pH 4.0-TTR in 1xLDS sample buffer (Life Technologies) was loaded onto a 10% NuPAGE bis-tris gel and subjected to electrophoresis in MES buffer at a constant 90V for 105 minutes. After electrophoresis, the gel was either stained in Instant Blue (Expedeon) or transferred to nitrocellulose filters for Western blot analysis.

10

#### *i. Native-PAGE*

**[0324]** Electrophoresis on native Tris-glycine gels was carried out as follows. 0.1-1  $\mu$ g TTR or pH 4.0-TTR in 1xTris-glycine sample buffer (Life Technologies) was loaded onto a 10-20% Tris-glycine gel and subjected to electrophoresis in 1x Native Tris-glycine running buffer at a constant 120V for 105 minutes. After electrophoresis, the gel was either stained in Instant Blue (Expedeon) or transferred to nitrocellulose filters for Western blot analysis.

15

#### *j. Western Blot*

**[0325]** SDS- or Native- PAGE gels were blotted onto nitrocellulose filter paper (iBlot, P7 Program) and blocked with blocking buffer (Licor) for 30 minutes. The filters were then incubated in 0.5  $\mu$ g/ml primary antibody in blocking buffer for 1 hour at room temperature (or over-night at 4°C), followed by three, 10 minutes washes with 1xTBS. The filters were placed in IRDye 800CW-conjugated goat-anti-mouse secondary diluted 1:20,000 in block buffer. After incubating the filters in secondary antibody solution for 1 hour at room temperature, the filters were washed and imaged on an Odyssey CLx infrared imager (Licor).

25

#### *k. TTR Fiber Formation Assay Protocol*

**[0326]** A solution of 3.6  $\mu$ M (0.2mg/ml) TTR-Y78F in 50mM sodium acetate, pH 4.8 was incubated at 37°C for 72 hours in the presence of 1.4  $\mu$ M (0.2mg/ml) mis-TTR antibody or an isotype control. After incubation, a 5X molar excess of thioflavin-T was added to the mixture and allowed to bind for 30 minutes. Fluorometric measurements were measured at an emissions wavelength of 480nm with an excitation wavelength set at 440nm. The 0% inhibition was set as the fluorescence intensity in the presence of an isotype control antibody (83 a.u.) and the 100% inhibition point was set as the fluorescence in the absence of TTR-Y78F protein (38 a.u.).

30

35

**[0327]** A stock solution of TTR-V122I (approximately 5 mg/mL) was diluted into buffer with final concentrations of 0.2 mg/mL TTR, 30 mM sodium acetate, 5 mM sodium phosphate, 100 mM KCl, 1 mM EDTA, 0.02% sodium azide, and varying concentrations of monoclonal antibody at pH 7.2. This sample was dialyzed against the same buffer at pH 4.5 for 3h at room temperature. Samples were then incubated at 37°C for 72 hours. After incubation a 5X molar excess of thioflavin-T was added to was added to the mixture and allowed to bind for 30 minutes. Thioflavin-T fluorescence was monitored using a Photon Technology International C60 spectrofluorimeter. The photomultiplier gain was varied and excitation and emission slit widths set to 2-4 nm to maximize signal to noise. Fluorescence measurements were made using 430 nm and 480 nm as excitation and emission wavelengths, respectively.

40

45

#### *1. Cardiac Tissue Samples*

**[0328]** Fresh frozen and paraffin-processed blocks of cardiac tissue with confirmed diagnoses of ATTR mutations were obtained from Dr. Merrill Benson at Indiana University. Samples included eight fresh frozen samples and six FFPE samples and each sample was diagnosed with either ATTR or some other cardiac amyloidosis. The diagnosis of the tissue was further confirmed at Prothena via IHC staining with antibodies to kappa and lambda light chains and amyloid A prior to characterization with the TTR antibodies.

50

#### *m. Immunohistochemistry*

**[0329]** Immunohistochemistry was performed on lightly paraformaldehyde-fixed, 10  $\mu$ m slide-mounted cryosections and on 5  $\mu$ m paraffin sections. The immunoperoxidase method was the principal detection system, which was performed on the Leica Bond Rx (Leica Biosystems, Buffalo Grove, IL) using the Bond Polymer Refine Detection Kit (DS980, Leica

55

Biosystems). The primary antibodies were incubated for one hour (according to concentrations in Table 2) followed by incubation with anti-mouse and anti-rabbit polymeric HRP-linker antibody conjugates. The staining was visualized with a DAB chromogen, which produced a brown deposit. The slides were counterstained with hematoxylin, dehydrated in an ascending series of alcohols, cleared in xylenes, and coverslipped with CytoSeal 60 (Richard Allen Scientific; Kalamazoo, MI). Negative control consisted of performing the entire immunohistochemical procedure on adjacent sections with a non-immune IgG isotype control or an omission of the primary antibody.

*n. Demonstration of Amyloid: Congo Red and Thioflavin T Staining*

**[0330]** Congo red stain was performed to demonstrate TTR amyloid in the tissue using a kit from American MasterTech (Lodi, California). The staining was performed according to the manufacturer's procedure. Slides were stained in the Congo Red solution for 1 hour followed by differentiation in 1% sodium hydroxide for approximately 15 seconds. The slides were then rinsed in running water, dehydrated through an alcohol series of increasing concentrations, and cleared through three changes of xylenes, and coverslipped with CytoSeal 60.

**[0331]** A modified Thioflavin T staining protocol (Schmidt et al 1995.) was employed to determine the presence of TTR amyloid in the tissue. Briefly, slides were counterstained with a Mayers hematoxylin, rinsed in running water and stained with a filtered solution of 0.015% Thioflavin T (T3516-25G; Sigma-Aldrich, St. Louis, MO) in 50% ethanol for ten minutes. The slides were then rinsed in running water and differentiated in 1% (v/v) acetic acid for 10 minutes and rinsed three times in water. The slides were allowed to air dry before being coverslipped with ProLong Gold (Life Technologies).

*o. Image Analysis*

**[0332]** Slides were imaged with either an Olympus BX61 microscope, Hamamatsu Nanoscope 2.0HT digital slide scanner, or a Leica SPE spectral confocal system. Images were collected and stored as TIFF files.

*p. Analysis of Human Plasma Samples by SDS-PAGE/Western*

**[0333]** Six plasma samples from patients confirmed for V30M ATTR (#21, #22, #23, #24, #25, #27) and 6 samples (Sample #11, #12, #15, #18, #19, #20) from normal subjects were obtained from M. Saraiva (Porto University, Portugal). Sample #C6 was a normal human serum sample obtained from a commercial source (BioreclamationIVT). These plasma samples were separated by SDS-PAGE and Western blotted with 9D5 or 5A1 as follows. A 1.4  $\mu$ l volume of plasma was diluted 1:8 into 1 $\times$ LDS sample buffer in the absence of reducing agent (Life Technologies). Samples were subjected to SDS-PAGE separation and Western blotted with 0.5  $\mu$ g/ml 9D5 or 5A1 as described previously.

*q. Analysis of Human Plasma Samples by MesoScale Discovery (MSD) Plate Assay*

**[0334]** 96-well MSD plates were coated with monoclonal antibody 6C1 at a concentration of 4  $\mu$ g/mL in PBS and incubated for 2 hours at room temperature with shaking, or overnight at 4°C. Plates were washed three times with 1 $\times$ TBST before being blocked with 3% MSD Blocker A solution, 150  $\mu$ L per well for 1 hour shaking. A 30  $\mu$ l per well volume of human plasma samples diluted 1:10 in a sample buffer comprised of 0.6% globulin-free bovine serum albumin, 1.5 mM monobasic sodium phosphate, 8mM dibasic sodium phosphate, 145 mM sodium chloride, 0.05% Triton X-405, and 0.05% thimerosal was added to the blocked MSD plates for 1 hour. Plates were washed 3 times with 1 $\times$ TBST. A 50  $\mu$ l per well volume of 1  $\mu$ g/ml sulfo-tagged detection antibody (either 8C3 total TTR antibody of the Dako polyclonal antibody) in sample buffer was added for 1 hr. at room temperature with shaking. Plates were washed three times with 1 $\times$ TBST followed by the addition of 150  $\mu$ l per well 1X Read Buffer T solution (Meso Scale Discovery). Plates were then read in the MSD Sector imager.

*r. Generation of an MSD Standard Curve*

**[0335]** In order to quantitate the amount of non-native, 6C1-reactive TTR protein present in human plasma samples, a MSD standard curve was generated using recombinant TTR-F87M/L110M as a 6C1-reactive TTR standard. This TTR variant contains two amino acid substitutions that prevent tetramer formation and keeps the protein in the monomer state (Jiang et al. (2001) Biochemistry 40, 11442-11452). As such, this TTR variant is recognized by all mis-TTR mAbs and is therefore well-suited for use as a reference standard in the MSD assay.

**[0336]** To generate the standard curve, 96-well MSD plates were coated with mis-TTR antibody 6C1 at a concentration of 4  $\mu$ g/mL in PBS and incubated for 2 hours at room temperature with shaking, or overnight at 4°C. Plates were washed three times with 1 $\times$ TBST before being blocked with 3% MSD Blocker A solution, 150  $\mu$ L per well for 1 hour shaking. The blocked plates were then treated for 1 hour with 50  $\mu$ l per well of 25  $\mu$ g/mL TTR-F87M/L110M serially diluted 1:5

with the last dilution being a buffer blank. Plates were washed 3 times with 1×TBST before the addition of a 50 μl per well volume of 1 μg/ml SulfoTag detection antibody (8C3-SulfoTag or Dako pAb-SulfoTag) for 1 hour at room temperature with shaking. Both 8C3 mAb and the Dako antibody were coupled to the SulfoTag and could be used at the detection antibody since they bound to total TTR and were not conformation specific.

5 [0337] After treatment with the detection antibody, plates were washed three times with a 150 μl per well volume of 1×TBST, followed by the addition of 150 μl per well 1× Read Buffer T (MSD). Plates were read in the MSD Sector imager and a resulting TTR F87M/L110M calibration curve was generated.

#### 10 s. Phagocytosis Assay

[0338] A 1-mg/mL sample of TTR-F87M/L110M, native TTR or low-pH aggregated TTR-V30M was amine coupled with pHrodo dye for 15 min at 37°C with a protein:dye ratio of ~15:4 according to the manufacturer's specifications (Thermo Scientific). Excess pHrodo label was removed by diafiltration in a spin concentrator with a 10K molecular weight cutoff (Pierce Thermo) and the pHrodo-TTR was resuspended in 1× PBS.

15 [0339] THP-1 human monocytes were cultured in cell culture media (RMPI, 10% low IgG serum, pen/strep). A 20-μg/mL aliquot of pHrodo-labeled TTR was separately pre-incubated with 40 μg/mL antibody at 37°C in cell culture media for 30 min prior to the addition of 5E+04 THP-1 cells in a 1:1 volumetric ratio. After tissue culture incubation (3 h), cells were washed with cell culture media three times, incubated in media for 10 min, then washed twice with and resuspended in FACS buffer (1% FBS in PBS). Red pHrodo fluorescence intensity was detected using Texas Red channel filters. Epifluorescence microscopy was carried out in a similar fashion. After FACS analysis, the remaining cells were transferred to glass chamber slides and imaged by inverted microscopy. Mean fluorescence intensities were automatically calculated by averaging the relative fluorescence intensities of each individual cell.

#### 25 Example 12. Antibody-dependent uptake of TTR by THP-1 cells

25 [0340] TTR-F87M/L110M was covalently labeled with the pH-sensitive fluorescent dye pHrodo. The pHrodo tag has minimal fluorescence under physiological pH, but fluorescence is enhanced upon engulfment into the low pH environments of endocytic vesicles and thus marks cellular uptake of tagged particles. THP-1 monocytes were added to pHrodo-tagged TTR (native or non-native, TTR-F87M/L110M) after treatment with either 14G8 or the isotype control antibody. Low levels of fluorescence were observed with either antibody incubated with native TTR, and after incubation of the control antibody with non-native TTR. Fluorescence was increased, however, after 14G8 incubation with non-native TTR (Figure 7A), suggesting that non-native TTR is not efficiently phagocytized under basal conditions, however the addition of mis-TTR antibodies specifically elicits phagocytosis of non-native TTR.

30 [0341] Dose-dependent phagocytosis of pHrodo-labeled, large aggregated fibrillar particles of TTR was also demonstrated for each of the mis-TTR mAbs (Figure 7B). Maximum antibody-dependent uptake was variable for each mis-TTR mAb (6C1>9D5~ 14G8> 5A1), reaching a plateau at mAb concentrations between 5-10 μg/mL. Variable antibody potencies may reflect isotype differences and associated changes in effector function among the four mis-TTR mAbs. Controls, including untreated cells or those treated with an IgG1 isotype control, did not demonstrate detectable or enhanced fluorescence, respectively.

#### 40 Example 13. Evaluation of mis-TTR Antibodies in Transgenic Mouse Model

[0342] *In vivo* studies are conducted in a humanized transgenic mouse model V30M hTTR (Inoue et al., (2008) Specific pathogen free conditions prevent transthyretin amyloidosis in mouse models. Transgenic Research 17:817-826) to assess the efficacy of anti-TTR antibodies in the binding and removal of aggregated hTTR.

45 [0343] Transgenic mice are bred using standard procedures and their circulating hTTR levels are assessed by ELISA. Mice with a serum level of 200-400 μg/ml of hTTR are used for subsequent efficacy studies. The first set of studies examine the natural deposition of hTTR in transgenic mice. Detection of hTTR deposits begins at 12 months of age and is repeated every 3-6 months thereafter. Once an acceptable level of aggregates is seen in transgenic mice, efficacy studies are initiated. Animals are divided into three treatment groups (n=10/ group) and treated weekly for four weeks with an IP dose of vehicle, control antibody (isotype control, 10 mpk) or an anti-hTTR antibody (10 mpk). One week after the last treatment the mice are euthanized, tissues collected and processed, and then stained to assess the number and size of remaining TTR deposits. Quantitative methods and statistics are employed to determine the degree of clearance seen among groups.

50 [0344] In an alternative approach, hTTR aggregates are prepared *in vitro* and then injected into the kidney of transgenic mice to seed the deposition of new aggregates. Applicant has determined that the injection of these preparations can expedite the deposition of new aggregates in a predictable manner. Based on these findings, animals are sedated, the left kidney exposed and pre-aggregated hTTR material injected into the cortex of the kidney. After a suitable recovery

period, mice are divided into three treatment groups (n=10/ group) and treated weekly for four-eight weeks with an IP dose of vehicle, control antibody (isotype control, 10 mpk) or an anti-hTTR antibody (10 mpk). One week after the last treatment the mice are euthanized, the kidneys collected and processed, and then stained to assess the number and size of TTR deposits. Quantitative methods and statistics are employed to determine the degree of change seen among groups.

**Example 14. Evaluation of mis-TTR Antibodies in a Matrigel Implant Model**

**[0345]** Applicant has determined that pre-aggregated hTTR can be suspended in Matrigel (BD Bioscience, Cat #354263), allowed to solidify and then placed subcutaneously in mice. At four weeks post implantation, the Matrigel implant maintained its structure and the aggregated hTTR was still present within the implant. Moreover, the implant was well tolerated by the mice and anti-hTTR antibodies were able to penetrate and bind to the aggregates suspended in the Matrigel. Based on these findings, an antibody efficacy study is conducted. Animals are sedated and an implant containing pre-aggregated hTTR suspended in Matrigel placed subcutaneously in mice. After a suitable recovery period, mice are divided into three treatment groups (n=10/ group) and treated weekly, for two-four weeks with an IP dose of vehicle, control antibody (isotype control, 10 mpk) or an anti-hTTR antibody (10 mpk). After the last treatment, the mice are euthanized, the skin containing the implant collected and processed, and then the amount of TTR deposits remaining assessed using histological and/ or biochemical methods. Quantitative analysis and statistics are employed to determine the degree of clearance seen among groups.

**SEQUENCE LISTING**

**[0346]**

<110> LIU, Yue NIJJAR, Tarlochan CHAKRABARTTY, Avijit

<120> ANTI-TRANSTHYRETIN ANTIBODIES

<130> 057450/473373

<150> 62/109, 002 <151> 2015-01-28

<150> 62/266,556 <151> 2015-12-11

<160> 118

<170> FastSEQ for Windows Version 4.0

<210> 1

<211> 124

<212> PRT

<213> Artificial Sequence

<220>

<223> m9D5VH

<400> 1

Glu Val Lys Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
 1 5 10 15  
 Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
 20 25 30  
 Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val  
 35 40 45  
 Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 Lys Gly Arg Phe Thr Phe Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 65 70 75 80  
 Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
 85 90 95  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 100 105 110  
 Val Trp Gly Thr Gly Thr Thr Val Thr Val Ser Ser  
 115 120

<210> 2  
 <211> 123  
 <212> PRT  
 20 <213> Artificial Sequence  
  
 <220>  
 <223> 1SEQ\_H  
  
 25 <400> 2

Glu Val Lys Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
 1 5 10 15  
 Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Thr Tyr  
 20 25 30  
 Thr Met Ser Trp Ala Arg Gln Thr Pro Glu Lys Lys Leu Glu Trp Val  
 35 40 45  
  
 35 Ala Tyr Ile Ser Lys Gly Gly Gly Ser Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 65 70 75 80  
 Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Leu Tyr Tyr Cys  
 85 90 95  
 Ala Arg Gly Ala Met Phe Gly Asn Asp Phe Lys Tyr Pro Met Asp Arg  
 100 105 110  
 Trp Gly Gln Gly Thr Ser Val Thr Val Ser Ser  
 115 120

<210> 3  
 <211> 124  
 <212> PRT  
 50 <213> Artificial Sequence  
  
 <220>  
 <223> BAC02114  
  
 55 <400> 3

5           Gln Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
           1                   5                   10                   15  
           Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
           20                   25                   30  
 10       Glu Met Asn Trp Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Trp Val  
           35                   40                   45  
           Ser Tyr Ile Thr Ser Ser Gly Ser Thr Ile Tyr Tyr Ala Asp Ser Val  
           50                   55                   60  
           Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Ser Leu Tyr  
           65                   70                   75                   80  
 15       Leu Gln Met Asn Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys  
           85                   90                   95  
           Ala Arg Gly Gly Gln Gly Ser Arg Tyr Tyr Tyr Tyr Tyr Gly Met Asp  
           100                   105                   110  
 20       Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser  
           115                   120

<210> 4  
 <211> 120  
 <212> PRT  
 20 <213> Artificial Sequence  
  
 <220>  
 <223> AAX82494  
  
 25 <400> 4

30           Gln Val Gln Leu Gln Glu Ser Gly Gly Gly Leu Val Lys Pro Gly Gly  
           1                   5                   10                   15  
           Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
           20                   25                   30  
           Gly Met Ser Trp Val Arg Gln Thr Pro Asp Lys Arg Leu Glu Trp Val  
           35                   40                   45  
           Ala Thr Ile Ser Ser Gly Gly Ser Tyr Thr Tyr Tyr Pro Asp Ser Val  
           50                   55                   60  
 35       Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
           65                   70                   75                   80  
           Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
           85                   90                   95  
 40       Ala Arg Leu Tyr Tyr Gly Tyr Arg Tyr Tyr Phe Asp Tyr Trp Gly Gln  
           100                   105                   110  
           Gly Thr Met Val Thr Val Ser Ser

115 120  
  
 45 <210> 5  
      <211> 124  
      <212> PRT  
      <213> Artificial Sequence  
  
 50 <220>  
      <223> Hu9D5VHv1  
  
      <400> 5

55

5           Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
           1                   5                   10                   15  
           Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
                   20                   25                   30  
 10       Thr Met Ser Trp Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Leu Val  
                   35                   40                   45  
           Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
                   50                   55                   60  
           Lys Gly Arg Phe Thr Phe Ser Arg Asp Asn Ala Lys Asn Ser Leu Tyr  
                   65                   70                   75                   80  
 15       Leu Gln Ser Asn Ser Leu Lys Ala Glu Asp Thr Ala Val Tyr Tyr Cys  
                   85                   90                   95  
           Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
                   100                   105                   110  
 20       Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser  
                   115                   120

<210> 6  
 <211> 124  
 <212> PRT  
 20 <213> Artificial Sequence  
  
 <220>  
 <223> Hu9D5VHv2  
  
 25 <400> 6

30           Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
           1                   5                   10                   15  
           Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
                   20                   25                   30  
           Thr Met Ser Trp Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Leu Val  
                   35                   40                   45  
           Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
                   50                   55                   60  
 35       Lys Gly Arg Phe Thr Phe Ser Arg Asp Asn Ala Lys Asn Ser Leu Tyr  
                   65                   70                   75                   80  
           Leu Gln Ser Asn Leu Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys  
                   85                   90                   95  
 40       Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
                   100                   105                   110  
           Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser  
                   115                   120

<210> 7  
 <211> 124  
 45 <212> PRT  
 <213> Artificial Sequence  
  
 <220>  
 <223> Hu9D5VHv2b  
  
 50 <400> 7

55

5           Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
             1                           5                           10                           15  
           Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
                   20                           25                           30  
 10       Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val  
                   35                           40                           45  
           Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
                   50                           55                           60  
           Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
                   65                           70                           75                           80  
 15       Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
                   85                           90                           95  
           Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
                   100                           105                           110  
 20       Val Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser  
                   115                           120

<210> 8  
 <211> 124  
 <212> PRT  
 20 <213> Artificial Sequence  
  
 <220>  
 <223> Hu9D5VHv3  
  
 25 <400> 8

30           Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
             1                           5                           10                           15  
           Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
                   20                           25                           30  
           Thr Met Ser Trp Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Trp Val  
                   35                           40                           45  
           Ser Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
                   50                           55                           60  
 35       Lys Gly Arg Phe Thr Phe Ser Arg Asp Asn Ala Lys Asn Ser Leu Tyr  
                   65                           70                           75                           80  
           Leu Gln Ser Asn Leu Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys  
                   85                           90                           95  
           Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
                   100                           105                           110  
 40       Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser  
                   115                           120

<210> 9  
 <211> 124  
 45 <212> PRT  
       <213> Artificial Sequence  
  
       <220>  
       <223> Hu9D5VHv3b  
  
 50 <400> 9

55           Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
             1                           5                           10                           15  
           Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
                   20                           25                           30  
           Thr Met Ser Trp Val Arg Gln Thr Pro Gly Lys Arg Leu Glu Leu Val  
                   35                           40                           45

Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 65 70 75 80  
 5 Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
 85 90 95  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 100 105 110  
 10 Val Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser  
 115 120

<210> 10  
 <211> 124  
 <212> PRT  
 15 <213> Artificial Sequence

<220>  
 <223> Hu9D5VHv4

20 <400> 10

Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
 1 5 10 15  
 Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
 20 25 30  
 25 Thr Met Ser Trp Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Trp Val  
 35 40 45  
 Ser Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 30 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Ser Leu Tyr  
 65 70 75 80  
 Leu Gln Ser Asn Leu Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys  
 85 90 95  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 100 105 110  
 35 Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser  
 115 120

<210> 11  
 <211> 124  
 40 <212> PRT  
 <213> Artificial Sequence

<220>  
 <223> Hu9D5VHv4b

45 <400> 11

50

55

```

    Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly
    1          5          10          15
    Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr
    20          25          30
5   Thr Met Ser Trp Val Arg Gln Ala Pro Gly Lys Arg Leu Glu Leu Val
    35          40          45
    Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val
    50          55          60
10  Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr
    65          70          75          80
    Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys
    85          90          95
    Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp
    100         105         110
15  Val Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser

```

115 120

```

<210> 12
20 <211> 124
    <212> PRT
    <213> Artificial Sequence

<220>
25 <223> Hu9D5VHv5

<400> 12

```

```

    Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly
    30 1          5          10          15
    Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr
    20          25          30
    Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val
    35 35          40          45
    Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val
    50          55          60
    Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr
    65          70          75          80
    Leu Gln Met Asn Leu Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys
    85          90          95
40  Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp
    100         105         110
    Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser
    115          120

```

```

45 <210> 13
    <211> 5
    <212> PRT
    <213> Artificial Sequence

50 <220>
    <223> 9D5 CDRH1 - Kabat

<400> 13

```

```

55 Ser Tyr Thr Met Ser
    1          5

```

5  
 <210> 14  
 <211> 17  
 <212> PRT  
 <213> Artificial Sequence

<220>  
 <223> 9D5 CDRH2 - Kabat

10  
 <400> 14

Glu	Ile	Ser	Asn	Ser	Gly	Asp	Thr	Thr	Tyr	Tyr	Pro	Asp	Thr	Val	Lys
1				5					10					15	
Gly															

15  
 <210> 15  
 <211> 15  
 <212> PRT  
 <213> Artificial Sequence

<220>  
 <223> 9D5 CDRH3 - Kabat

20  
 <400> 15

His	Tyr	Tyr	Tyr	Gly	Gly	Gly	Tyr	Gly	Gly	Trp	Phe	Phe	Asp	Val
1				5					10					15

30  
 <210> 16  
 <211> 112  
 <212> PRT  
 <213> Artificial Sequence

<220>  
 <223> m9D5VL

35  
 <400> 16

Asp	Ile	Val	Met	Thr	Gln	Ala	Ala	Pro	Ser	Val	Pro	Val	Thr	Pro	Gly
1				5					10					15	
Glu	Ser	Val	Ser	Ile	Ser	Cys	Arg	Ser	Ser	Lys	Ser	Leu	Leu	His	Ser
			20					25					30		
Asn	Gly	Asn	Thr	Tyr	Leu	Tyr	Trp	Phe	Leu	Gln	Arg	Pro	Gly	Gln	Ser
			35				40					45			
Pro	Gln	Leu	Leu	Ile	Tyr	Arg	Val	Ser	Asn	Leu	Ala	Ser	Gly	Val	Pro
			50			55					60				
Asp	Arg	Phe	Ser	Gly	Ser	Gly	Ser	Gly	Thr	Ala	Phe	Thr	Leu	Arg	Ile
					70					75				80	
Ser	Arg	Val	Glu	Ala	Glu	Asp	Val	Gly	Val	Tyr	Tyr	Cys	Met	Gln	His
				85					90					95	
Leu	Glu	Tyr	Pro	Leu	Thr	Phe	Gly	Ala	Gly	Thr	Lys	Leu	Glu	Leu	Lys
			100					105					110		

55  
 <210> 17  
 <211> 112  
 <212> PRT  
 <213> Artificial Sequence

<220>  
<223> 1MJU\_L

<400> 17

5

10

15

20

```

Asp Ile Val Met Thr Gln Ala Ala Pro Ser Val Pro Val Thr Pro Gly
 1           5           10           15
Glu Ser Val Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
           20           25           30
Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg Pro Gly Gln Ser
           35           40           45
Pro Gln Leu Leu Ile Tyr Arg Met Ser Asn Leu Ala Ser Gly Val Pro
           50           55           60
Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Arg Ile
65           70           75           80
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Leu Gln His
           85           90           95
Leu Glu Tyr Pro Phe Thr Phe Gly Ala Gly Thr Lys Leu Glu Leu Lys
           100          105          110
    
```

<210> 18  
<211> 112  
<212> PRT  
<213> Artificial Sequence

25

<220>  
<223> ABC66952

<400> 18

30

35

40

```

Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
 1           5           10           15
Glu Pro Ala Ser Ile Ser Cys Arg Ser Asn Gln Ser Leu Leu Tyr Ser
           20           25           30
Asn Gly Tyr Asn Tyr Leu Asp Trp Tyr Leu Gln Lys Pro Gly Gln Ser
           35           40           45
Pro Gln Leu Leu Ile Tyr Ser Gly Ser Asn Arg Ala Ser Gly Val Pro
           50           55           60
Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
65           70           75           80
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln Ala
           85           90           95
Leu Gln Ser Pro Tyr Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
           100          105          110
    
```

45

<210> 19  
<211> 112  
<212> PRT  
<213> Artificial Sequence

50

<220>  
<223> Hu9D5VLv1

<400> 19

55

```

    Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
    1           5           10           15
    Glu Pro Ala Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
    20           25           30
5   Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Lys Pro Gly Gln Ser
    35           40           45
    Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
    50           55           60
10  Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
    65           70           75           80
    Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
    85           90           95
    Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
    100          105          110

```

```

15  <210> 20
    <211> 112
    <212> PRT
    <213> Artificial Sequence
20
    <220>
    <223> Hu9D5VLv2
25
    <400> 20

```

```

    Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
    1           5           10           15
    Glu Pro Ala Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
    20           25           30
30  Asn Gly Asn Thr Tyr Leu Tyr Trp Tyr Leu Gln Lys Pro Gly Gln Ser
    35           40           45
    Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
    50           55           60
35  Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
    65           70           75           80
    Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
    85           90           95
40  Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
    100          105          110

```

```

45  <210> 21
    <211> 112
    <212> PRT
    <213> Artificial Sequence
    <220>
    <223> Hu9D5VLv3
50
    <400> 21

```

55

Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly  
 1 5 10 15  
 Glu Pro Ala Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser  
 20 25 30  
 5 Asn Gly Asn Thr Tyr Leu Tyr Trp Tyr Leu Gln Lys Pro Gly Gln Ser  
 35 40 45  
 Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro  
 50 55 60  
 Ser Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile  
 65 70 75 80  
 Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His  
 85 90 95  
 Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys  
 100 105 110

15  
 <210> 22  
 <211> 112  
 <212> PRT  
 <213> Artificial Sequence

20  
 <220>  
 <223> Hu9D5VLv4

25  
 <400> 22

Asp Ile Val Met Thr Gln Ser Ala Pro Ser Leu Pro Val Thr Pro Gly  
 1 5 10 15  
 Glu Pro Val Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser  
 20 25 30  
 30 Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg Pro Gly Gln Ser  
 35 40 45  
 Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro  
 50 55 60  
 Ser Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Arg Ile  
 65 70 75 80  
 Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His  
 85 90 95  
 Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys  
 100 105 110

40  
 <210> 23  
 <211> 112  
 <212> PRT  
 <213> Artificial Sequence

45  
 <220>  
 <223> Hu9D5VLv5

50  
 <400> 23

55

```

    Asp Ile Val Met Thr Gln Ser Ala Pro Ser Leu Pro Val Thr Pro Gly
    1           5           10           15
    Glu Ser Val Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
    20           25           30
5   Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg Pro Gly Gln Ser
    35           40           45
    Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
    50           55           60
    Ser Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Arg Ile
    65           70           75           80
    Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
    85           90           95
    Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
    100          105          110

```

```

15  <210> 24
    <211> 16
    <212> PRT
    <213> Artificial Sequence

```

```

20  <220>
    <223> 9D5 CDRL1 - Kabat

```

```

25  <400> 24

```

```

    Arg Ser Ser Lys Ser Leu Leu His Ser Asn Gly Asn Thr Tyr Leu Tyr
    1           5           10           15

```

```

30  <210> 25
    <211> 7
    <212> PRT
    <213> Artificial Sequence

```

```

35  <220>
    <223> 9D5 CDRL2 - Kabat

```

```

    <400> 25

```

```

    Arg Val Ser Asn Leu Ala Ser
    1           5

```

```

40  <210> 26
    <211> 9
    <212> PRT
    <213> Artificial Sequence

```

```

45  <220>
    <223> 9D5 CDRL3 - Kabat

```

```

50  <400> 26

```

```

    Met Gln His Leu Glu Tyr Pro Leu Thr
    1           5

```

```

55  <210> 27
    <211> 454
    <212> PRT
    <213> Artificial Sequence

```

<220>

<223> Hu9D5 H1 heavy chain

<400> 27

5

10

15

20

25

30

35

40

45

50

55

	Glu	Val	Gln	Leu	Val	Glu	Ser	Gly	Gly	Gly	Leu	Val	Gln	Pro	Gly	Gly
	1				5					10					15	
	Ser	Leu	Arg	Leu	Ser	Cys	Ala	Ala	Ser	Gly	Phe	Thr	Phe	Ser	Ser	Tyr
5				20					25					30		
	Thr	Met	Ser	Trp	Val	Arg	Gln	Ala	Pro	Gly	Lys	Gly	Leu	Glu	Leu	Val
		35						40					45			
	Ala	Glu	Ile	Ser	Asn	Ser	Gly	Asp	Thr	Thr	Tyr	Tyr	Pro	Asp	Thr	Val
		50					55					60				
10	Lys	Gly	Arg	Phe	Thr	Phe	Ser	Arg	Asp	Asn	Ala	Lys	Asn	Ser	Leu	Tyr
	65					70					75					80
	Leu	Gln	Ser	Asn	Ser	Leu	Lys	Ala	Glu	Asp	Thr	Ala	Val	Tyr	Tyr	Cys
				85						90					95	
	Ala	Arg	His	Tyr	Tyr	Tyr	Gly	Gly	Gly	Tyr	Gly	Gly	Trp	Phe	Phe	Asp
			100						105					110		
15	Val	Trp	Gly	Gln	Gly	Thr	Thr	Val	Thr	Val	Ser	Ser	Ala	Ser	Thr	Lys
			115					120					125			
	Gly	Pro	Ser	Val	Phe	Pro	Leu	Ala	Pro	Ser	Ser	Lys	Ser	Thr	Ser	Gly
		130					135					140				
	Gly	Thr	Ala	Ala	Leu	Gly	Cys	Leu	Val	Lys	Asp	Tyr	Phe	Pro	Glu	Pro
20		145				150					155					160
	Val	Thr	Val	Ser	Trp	Asn	Ser	Gly	Ala	Leu	Thr	Ser	Gly	Val	His	Thr
				165						170					175	
	Phe	Pro	Ala	Val	Leu	Gln	Ser	Ser	Gly	Leu	Tyr	Ser	Leu	Ser	Ser	Val
			180						185				190			
25	Val	Thr	Val	Pro	Ser	Ser	Ser	Leu	Gly	Thr	Gln	Thr	Tyr	Ile	Cys	Asn
			195					200					205			
	Val	Asn	His	Lys	Pro	Ser	Asn	Thr	Lys	Val	Asp	Lys	Arg	Val	Glu	Pro
		210					215					220				
	Lys	Ser	Cys	Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala	Pro	Glu
		225				230					235					240
30	Leu	Leu	Gly	Gly	Pro	Ser	Val	Phe	Leu	Phe	Pro	Pro	Lys	Pro	Lys	Asp
				245						250					255	
	Thr	Leu	Met	Ile	Ser	Arg	Thr	Pro	Glu	Val	Thr	Cys	Val	Val	Val	Asp
			260						265						270	
	Val	Ser	His	Glu	Asp	Pro	Glu	Val	Lys	Phe	Asn	Trp	Tyr	Val	Asp	Gly
			275				280						285			
35	Val	Gln	Val	His	Asn	Ala	Lys	Thr	Lys	Pro	Arg	Glu	Glu	Gln	Tyr	Asn
		290					295					300				
	Ser	Thr	Tyr	Arg	Val	Val	Ser	Val	Leu	Thr	Val	Leu	His	Gln	Asp	Trp
				305			310				315					320
	Leu	Asn	Gly	Lys	Glu	Tyr	Lys	Cys	Lys	Val	Ser	Asn	Lys	Ala	Leu	Pro
40				325							330					335
	Ala	Pro	Ile	Glu	Lys	Thr	Ile	Ser	Lys	Ala	Lys	Gly	Gln	Pro	Arg	Glu
			340						345						350	
	Pro	Gln	Val	Tyr	Thr	Leu	Pro	Pro	Ser	Arg	Glu	Glu	Met	Thr	Lys	Asn
			355						360					365		
45	Gln	Val	Ser	Leu	Thr	Cys	Leu	Val	Lys	Gly	Phe	Tyr	Pro	Ser	Asp	Ile
		370					375					380				
	Ala	Val	Glu	Trp	Glu	Ser	Asn	Gly	Gln	Pro	Glu	Asn	Asn	Tyr	Lys	Thr
		385				390					395					400
	Thr	Pro	Pro	Val	Leu	Asp	Ser	Asp	Gly	Ser	Phe	Phe	Leu	Tyr	Ser	Lys
				405						410					415	
50	Leu	Thr	Val	Asp	Lys	Ser	Arg	Trp	Gln	Gln	Gly	Asn	Val	Phe	Ser	Cys
			420						425						430	
	Ser	Val	Met	His	Glu	Ala	Leu	His	Asn	His	Tyr	Thr	Gln	Lys	Ser	Leu
			435					440						445		
55	Ser	Leu	Ser	Pro	Gly	Lys										
		450														

<210> 28  
<211> 454

&lt;212&gt; PRT

&lt;213&gt; Artificial Sequence

&lt;220&gt;

5 &lt;223&gt; Hu9D5 H2 heavy chain

&lt;400&gt; 28

```

10      Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly
        1           5           10           15
        Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr
        20
        Thr Met Ser Trp Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Leu Val
        35           40           45
15      Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val
        50           55           60
        Lys Gly Arg Phe Thr Phe Ser Arg Asp Asn Ala Lys Asn Ser Leu Tyr
        65           70           75           80
        Leu Gln Ser Asn Leu Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys
        85           90           95
20      Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp
        100          105          110
        Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser Ala Ser Thr Lys
        115          120          125
25      Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys Ser Thr Ser Gly
        130          135          140
        Gly Thr Ala Ala Leu Gly Cys Leu Val Lys Asp Tyr Phe Pro Glu Pro
        145          150          155          160
        Val Thr Val Ser Trp Asn Ser Gly Ala Leu Thr Ser Gly Val His Thr
        165          170          175
30      Phe Pro Ala Val Leu Gln Ser Ser Gly Leu Tyr Ser Leu Ser Ser Val
        180          185          190
        Val Thr Val Pro Ser Ser Ser Leu Gly Thr Gln Thr Tyr Ile Cys Asn
        195          200          205
        Val Asn His Lys Pro Ser Asn Thr Lys Val Asp Lys Arg Val Glu Pro
        210          215          220
35      Lys Ser Cys Asp Lys Thr His Thr Cys Pro Pro Cys Pro Ala Pro Glu
        225          230          235          240
        Leu Leu Gly Gly Pro Ser Val Phe Leu Phe Pro Pro Lys Pro Lys Asp
        245          250          255          260
40      Thr Leu Met Ile Ser Arg Thr Pro Glu Val Thr Cys Val Val Val Asp
        260          265          270
        Val Ser His Glu Asp Pro Glu Val Lys Phe Asn Trp Tyr Val Asp Gly
        275          280          285
        Val Glu Val His Asn Ala Lys Thr Lys Pro Arg Glu Glu Gln Tyr Asn
        290          295          300
45      Ser Thr Tyr Arg Val Val Ser Val Leu Thr Val Leu His Gln Asp Trp
        305          310          315          320
        Leu Asn Gly Lys Glu Tyr Lys Cys Lys Val Ser Asn Lys Ala Leu Pro
        325          330          335
        Ala Pro Ile Glu Lys Thr Ile Ser Lys Ala Lys Gly Gln Pro Arg Glu
        340          345          350
50      Pro Gln Val Tyr Thr Leu Pro Pro Ser Arg Glu Glu Met Thr Lys Asn
        355          360          365
        Gln Val Ser Leu Thr Cys Leu Val Lys Gly Phe Tyr Pro Ser Asp Ile
        370          375          380
        Ala Val Glu Trp Glu Ser Asn Gly Gln Pro Glu Asn Asn Tyr Lys Thr
        385          390          395          400
55      Thr Pro Pro Val Leu Asp Ser Asp Gly Ser Phe Phe Leu Tyr Ser Lys
        405          410          415
        Leu Thr Val Asp Lys Ser Arg Trp Gln Gln Gly Asn Val Phe Ser Cys

```

```

                420                425                430
Ser Val Met His Gln Ala Leu His Asn His Tyr Thr Gln Lys Ser Leu
                435                440                445
Ser Leu Ser Pro Gly Lys
5

```

<210> 29

<211> 454

<212> PRT

10 <213> Artificial Sequence

<220>

<223> Hu9D5 H2b heavy chain

15 <400> 29

20

25

30

35

40

45

50

55

Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
 1 5 10  
 Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
 20 25 30  
 Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val  
 35 40 45  
 Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 65 70 75 80  
 Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
 85 90 95  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 100 105 110  
 Val Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser Ala Ser Thr Lys  
 115 120 125  
 Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys Ser Thr Ser Gly  
 130 135 140  
 Gly Thr Ala Ala Leu Gly Cys Leu Val Lys Asp Tyr Phe Pro Glu Pro  
 145 150 155 160  
 Val Thr Val Ser Trp Asn Ser Gly Ala Leu Thr Ser Gly Val His Thr  
 165 170 175  
 Phe Pro Ala Val Leu Gln Ser Ser Gly Leu Tyr Ser Leu Ser Ser Val  
 180 185 190  
 Val Thr Val Pro Ser Ser Ser Leu Gly Thr Gln Thr Tyr Ile Cys Asn  
 195 200 205  
 Val Asn His Lys Pro Ser Asn Thr Lys Val Asp Lys Arg Val Glu Pro  
 210 215 220  
 Lys Ser Cys Asp Lys Thr His Thr Cys Pro Pro Cys Pro Ala Pro Glu  
 225 230 235 240  
 Leu Leu Gly Gly Pro Ser Val Phe Leu Phe Pro Pro Lys Pro Lys Asp  
 245 250 255  
 Thr Leu Met Ile Ser Arg Thr Pro Glu Val Thr Cys Val Val Val Asp  
 260 265 270  
 Val Ser His Glu Asp Pro Glu Val Lys Phe Asn Trp Tyr Val Asp Gly  
 275 280 285  
 Val Glu Val His Asn Ala Lys Thr Lys Pro Arg Glu Glu Gln Tyr Asn  
 290 295 300  
 Ser Thr Tyr Arg Val Val Ser Val Leu Thr Val Leu His Gln Asp Trp  
 305 310 315 320  
 Leu Asn Gly Lys Glu Tyr Lys Cys Lys Val Ser Asn Lys Ala Leu Pro  
 325 330 335  
 Ala Pro Ile Glu Lys Thr Ile Ser Lys Ala Lys Gly Gln Pro Arg Glu  
 340 345 350  
 Pro Gln Val Tyr Thr Leu Pro Pro Ser Arg Glu Glu Met Thr Lys Asn  
 355 360 365  
 Gln Val Ser Leu Thr Cys Leu Val Lys Gly Phe Tyr Pro Ser Asp Ile  
 370 375 380  
 Ala Val Glu Trp Glu Ser Asn Gly Gln Pro Glu Asn Asn Tyr Lys Thr  
 385 390 395 400  
 Thr Pro Pro Val Leu Asp Ser Asp Gly Ser Phe Phe Leu Tyr Ser Lys  
 405 410 415  
 Leu Thr Val Asp Lys Ser Arg Trp Gln Gln Gly Asn Val Phe Ser Cys  
 420 425 430  
 Ser Val Met His Glu Ala Leu His Asn His Tyr Thr Gln Lys Ser Leu  
 435 440 445  
 Ser Leu Ser Pro Gly Lys  
 450

<210> 30  
 <211> 454  
 <212> PRT  
 <213> Artificial Sequence

5

<220>  
 <223> Hu9D5 H3 heavy chain

10

&lt;400&gt; 30

	Glu	Val	Gln	Leu	Val	Glu	Ser	Gly	Gly	Gly	Leu	Val	Gln	Pro	Gly	Gly
	1				5					10					15	
	Ser	Leu	Arg	Leu	Ser	Cys	Ala	Ala	Ser	Gly	Phe	Thr	Phe	Ser	Ser	Tyr
				20					25					30		
15	Thr	Met	Ser	Trp	Val	Arg	Gln	Ala	Pro	Gly	Lys	Gly	Leu	Glu	Trp	Val
			35					40					45			
	Ser	Glu	Ile	Ser	Asn	Ser	Gly	Asp	Thr	Thr	Tyr	Tyr	Pro	Asp	Thr	Val
	50					55						60				
	Lys	Gly	Arg	Phe	Thr	Phe	Ser	Arg	Asp	Asn	Ala	Lys	Asn	Ser	Leu	Tyr
20	65					70					75					80
	Leu	Gln	Ser	Asn	Leu	Leu	Arg	Ala	Glu	Asp	Thr	Ala	Val	Tyr	Tyr	Cys
				85						90					95	
	Ala	Arg	His	Tyr	Tyr	Tyr	Gly	Gly	Gly	Tyr	Gly	Gly	Trp	Phe	Phe	Asp
			100						105					110		
25	Val	Trp	Gly	Gln	Gly	Thr	Thr	Val	Thr	Val	Ser	Ser	Ala	Ser	Thr	Lys
			115					120						125		
	Gly	Pro	Ser	Val	Phe	Pro	Leu	Ala	Pro	Ser	Ser	Lys	Ser	Thr	Ser	Gly
	130						135						140			
	Gly	Thr	Ala	Ala	Leu	Gly	Cys	Leu	Val	Lys	Asp	Tyr	Phe	Pro	Glu	Pro
	145					150					155					160
30	Val	Thr	Val	Ser	Trp	Asn	Ser	Gly	Ala	Leu	Thr	Ser	Gly	Val	His	Thr
				165							170				175	
	Phe	Pro	Ala	Val	Leu	Gln	Ser	Ser	Gly	Leu	Tyr	Ser	Leu	Ser	Ser	Val
			180						185					190		
	Val	Thr	Val	Pro	Ser	Ser	Ser	Leu	Gly	Thr	Gln	Thr	Tyr	Ile	Cys	Asn
			195					200					205			
35	Val	Asn	His	Lys	Pro	Ser	Asn	Thr	Lys	Val	Asp	Lys	Arg	Val	Glu	Pro
	210						215					220				
	Lys	Ser	Cys	Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala	Pro	Glu
	225					230					235					240
40	Leu	Leu	Gly	Gly	Pro	Ser	Val	Phe	Leu	Phe	Pro	Pro	Lys	Pro	Lys	Asp
				245					250					255		
	Thr	Leu	Met	Ile	Ser	Arg	Thr	Pro	Glu	Val	Thr	Cys	Val	Val	Val	Asp
			260						265					270		
	Val	Ser	His	Glu	Asp	Pro	Glu	Val	Lys	Phe	Asn	Trp	Tyr	Val	Asp	Gly
			275				280						285			
45	Val	Glu	Val	His	Asn	Ala	Lys	Thr	Lys	Pro	Arg	Glu	Glu	Gln	Tyr	Asn
	290					295						300				
	Ser	Thr	Tyr	Arg	Val	Val	Ser	Val	Leu	Thr	Val	Leu	His	Gln	Asp	Trp
	305					310					315				320	
50	Leu	Asn	Gly	Lys	Glu	Tyr	Lys	Cys	Lys	Val	Ser	Asn	Lys	Ala	Leu	Pro

50

55

					325					330					335	
	Ala	Pro	Ile	Glu	Lys	Thr	Ile	Ser	Lys	Ala	Lys	Gly	Gln	Pro	Arg	Glu
				340					345					350		
5	Pro	Gln	Val	Tyr	Thr	Leu	Pro	Pro	Ser	Arg	Glu	Glu	Met	Thr	Lys	Asn
			355					360					365			
	Gln	Val	Ser	Leu	Thr	Cys	Leu	Val	Lys	Gly	Phe	Tyr	Pro	Ser	Asp	Ile
		370				375					380					
	Ala	Val	Glu	Trp	Glu	Ser	Asn	Gly	Gln	Pro	Glu	Asn	Asn	Tyr	Lys	Thr
10		385				390					395					400
	Thr	Pro	Pro	Val	Leu	Asp	Ser	Asp	Gly	Ser	Phe	Phe	Leu	Tyr	Ser	Lys
				405					410						415	
	Leu	Thr	Val	Asp	Lys	Ser	Arg	Trp	Gln	Gln	Gly	Asn	Val	Phe	Ser	Cys
			420						425					430		
	Ser	Val	Met	His	Glu	Ala	Leu	His	Asn	His	Tyr	Thr	Gln	Lys	Ser	Leu
15			435					440					445			
	Ser	Leu	Ser	Pro	Gly	Lys										
			450													

<210> 31  
 <211> 454  
 <212> PRT  
 <213> Artificial Sequence  
  
 <220>  
 <223> Hu9D5 H3b heavy chain  
  
 <400> 31

30

35

40

45

50

55

Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
 1 5 10  
 Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
 20 25 30  
 5 Thr Met Ser Trp Val Arg Gln Thr Pro Gly Lys Arg Leu Glu Leu Val  
 35 40 45  
 Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 10 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 65 70 75 80  
 Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
 85 90 95  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 100 105 110  
 15 Val Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser Ala Ser Thr Lys  
 115 120 125  
 Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys Ser Thr Ser Gly  
 130 135 140  
 Gly Thr Ala Ala Leu Gly Cys Leu Val Lys Asp Tyr Phe Pro Glu Pro  
 145 150 155 160  
 20 Val Thr Val Ser Trp Asn Ser Gly Ala Leu Thr Ser Gly Val His Thr  
 165 170 175  
 Phe Pro Ala Val Leu Gln Ser Ser Gly Leu Tyr Ser Leu Ser Ser Val  
 180 185 190  
 Val Thr Val Pro Ser Ser Ser Leu Gly Thr Gln Thr Tyr Ile Cys Asn  
 195 200 205  
 25 Val Asn His Lys Pro Ser Asn Thr Lys Val Asp Lys Arg Val Glu Pro  
 210 215 220  
 Lys Ser Cys Asp Lys Thr His Thr Cys Pro Pro Cys Pro Ala Pro Glu  
 225 230 235 240  
 30 Leu Leu Gly Gly Pro Ser Val Phe Leu Phe Pro Pro Lys Pro Lys Asp  
 245 250 255  
 Thr Leu Met Ile Ser Arg Thr Pro Glu Val Thr Cys Val Val Val Asp  
 260 265 270  
 Val Ser His Glu Asp Pro Glu Val Lys Phe Asn Trp Tyr Val Asp Gly  
 275 280 285  
 35 Val Glu Val His Asn Ala Lys Thr Lys Pro Arg Glu Glu Gln Tyr Asn  
 290 295 300  
 Ser Thr Tyr Arg Val Val Ser Val Leu Thr Val Leu His Gln Asp Trp  
 305 310 315 320  
 40 Leu Asn Gly Lys Glu Tyr Lys Cys Lys Val Ser Asn Lys Ala Leu Pro  
 325 330 335  
 Ala Pro Ile Glu Lys Thr Ile Ser Lys Ala Lys Gly Gln Pro Arg Glu  
 340 345 350  
 Pro Gln Val Tyr Thr Leu Pro Pro Ser Arg Glu Glu Met Thr Lys Asn  
 355 360 365  
 45 Gln Val Ser Leu Thr Cys Leu Val Lys Gly Phe Tyr Pro Ser Asp Ile  
 370 375 380  
 Ala Val Glu Trp Glu Ser Asn Gly Gln Pro Glu Asn Asn Tyr Lys Thr  
 385 390 395 400  
 50 Thr Pro Pro Val Leu Asp Ser Asp Gly Ser Phe Phe Leu Tyr Ser Lys  
 405 410 415  
 Leu Thr Val Asp Lys Ser Arg Trp Gln Gln Gly Asn Val Phe Ser Cys  
 420 425 430  
 Ser Val Met His Glu Ala Leu His Asn His Tyr Thr Gln Lys Ser Leu  
 435 440 445  
 55 Ser Leu Ser Pro Gly Lys  
 450

<211> 454  
 <212> PRT  
 <213> Artificial Sequence

5 <220>  
 <223> Hu9D5 H4 heavy chain

<400> 32

10           Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
               1                  5                  10                  15  
               Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
                           20                  25                  30  
 15           Thr Met Ser Trp Val Arg Gln Ala Pro Gly Lys Gly Leu Glu Trp Val  
                   35                  40                  45  
               Ser Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
                   50                  55                  60  
               Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Ser Leu Tyr  
                   65                  70                  75                  80  
 20           Leu Gln Ser Asn Leu Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys  
                   85                  90                  95  
               Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
                   100                  105                  110  
               Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser Ala Ser Thr Lys  
                   115                  120                  125  
 25           Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys Ser Thr Ser Gly  
                   130                  135                  140  
               Gly Thr Ala Ala Leu Gly Cys Leu Val Lys Asp Tyr Phe Pro Glu Pro  
                   145                  150                  155                  160  
               Val Thr Val Ser Trp Asn Ser Gly Ala Leu Thr Ser Gly Val His Thr  
                   165                  170                  175  
 30           Phe Pro Ala Val Leu Gln Ser Ser Gly Leu Tyr Ser Leu Ser Ser Val  
                   180                  185                  190  
               Val Thr Val Pro Ser Ser Ser Leu Gly Thr Gln Thr Tyr Ile Cys Asn  
                   195                  200                  205  
 35           Val Asn His Lys Pro Ser Asn Thr Lys Val Asp Lys Arg Val Glu Pro  
                   210                  215                  220  
               Lys Ser Cys Asp Lys Thr His Thr Cys Pro Pro Cys Pro Ala Pro Glu

40

45

50

55

	225				230				235				240			
	Leu	Leu	Gly	Gly	Pro	Ser	Val	Phe	Leu	Phe	Pro	Pro	Lys	Pro	Lys	Asp
					245				250				255			
5	Thr	Leu	Met	Ile	Ser	Arg	Thr	Pro	Glu	Val	Thr	Cys	Val	Val	Val	Asp
			260					265					270			
	Val	Ser	His	Glu	Asp	Pro	Glu	Val	Lys	Phe	Asn	Trp	Tyr	Val	Asp	Gly
			275					280					285			
	Val	Glu	Val	His	Asn	Ala	Lys	Thr	Lys	Pro	Arg	Glu	Glu	Gln	Tyr	Asn
			290					295					300			
10	Ser	Thr	Tyr	Arg	Val	Val	Ser	Val	Leu	Thr	Val	Leu	His	Gln	Asp	Trp
			305				310					315				320
	Leu	Asn	Gly	Lys	Glu	Tyr	Lys	Cys	Lys	Val	Ser	Asn	Lys	Ala	Leu	Pro
			325							330						335
	Ala	Pro	Ile	Glu	Lys	Thr	Ile	Ser	Lys	Ala	Lys	Gly	Gln	Pro	Arg	Glu
			340							345						350
15	Pro	Gln	Val	Tyr	Thr	Leu	Pro	Pro	Ser	Arg	Glu	Glu	Met	Thr	Lys	Asn
			355							360						365
	Gln	Val	Ser	Leu	Thr	Cys	Leu	Val	Lys	Gly	Phe	Tyr	Pro	Ser	Asp	Ile
			370							375						380
	Ala	Val	Glu	Trp	Glu	Ser	Asn	Gly	Gln	Pro	Glu	Asn	Asn	Tyr	Lys	Thr
			385				390					395				400
20	Thr	Pro	Pro	Val	Leu	Asp	Ser	Asp	Gly	Ser	Phe	Phe	Leu	Tyr	Ser	Lys
						405										415
	Leu	Thr	Val	Asp	Lys	Ser	Arg	Trp	Gln	Gln	Gly	Asn	Val	Phe	Ser	Cys
			420						425							430
25	Ser	Val	Met	His	Glu	Ala	Leu	His	Asn	His	Tyr	Thr	Gln	Lys	Ser	Leu
			435						440							445
	Ser	Leu	Ser	Pro	Gly	Lys										
			450													

30 <210> 33  
 <211> 454  
 <212> PRT  
 <213> Artificial Sequence

35 <220>  
 <223> Hu9D5 H4b heavy chain

<400> 33

40

45

50

55

	Glu	Val	Gln	Leu	Val	Glu	Ser	Gly	Gly	Gly	Leu	Val	Gln	Pro	Gly	Gly
	1				5					10					15	
	Ser	Leu	Lys	Leu	Ser	Cys	Ala	Ala	Ser	Gly	Phe	Thr	Phe	Ser	Ser	Tyr
5				20					25					30		
	Thr	Met	Ser	Trp	Val	Arg	Gln	Ala	Pro	Gly	Lys	Arg	Leu	Glu	Leu	Val
			35					40					45			
	Ala	Gln	Ile	Ser	Asn	Ser	Gly	Asp	Thr	Thr	Tyr	Tyr	Pro	Asp	Thr	Val
	50						55					60				
10	Lys	Gly	Arg	Phe	Thr	Ile	Ser	Arg	Asp	Asn	Ala	Lys	Asn	Thr	Leu	Tyr
	65					70					75					80
	Leu	Gln	Met	Ser	Ser	Leu	Lys	Ser	Glu	Asp	Thr	Ala	Met	Tyr	Tyr	Cys
						85				90					95	
	Ala	Arg	His	Tyr	Tyr	Tyr	Gly	Gly	Gly	Tyr	Gly	Gly	Trp	Phe	Phe	Asp
				100					105					110		
15	Val	Trp	Gly	Gln	Gly	Thr	Leu	Val	Thr	Val	Ser	Ser	Ala	Ser	Thr	Lys
			115					120					125			
	Gly	Pro	Ser	Val	Phe	Pro	Leu	Ala	Pro	Ser	Ser	Lys	Ser	Thr	Ser	Gly
	130						135					140				
	Gly	Thr	Ala	Ala	Leu	Gly	Cys	Leu	Val	Lys	Asp	Tyr	Phe	Pro	Glu	Pro
20	145					150					155					160
	Val	Thr	Val	Ser	Trp	Asn	Ser	Gly	Ala	Leu	Thr	Ser	Gly	Val	His	Thr
					165					170					175	
	Phe	Pro	Ala	Val	Leu	Gln	Ser	Ser	Gly	Leu	Tyr	Ser	Leu	Ser	Ser	Val
25				180					185					190		
	Val	Thr	Val	Pro	Ser	Ser	Ser	Leu	Gly	Thr	Gln	Thr	Tyr	Ile	Cys	Asn
			195					200					205			
	Val	Asn	His	Lys	Pro	Ser	Asn	Thr	Lys	Val	Asp	Lys	Arg	Val	Glu	Pro
	210						215					220				
30	Lys	Ser	Cys	Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala	Pro	Glu
	225					230					235					240
	Leu	Leu	Gly	Gly	Pro	Ser	Val	Phe	Leu	Phe	Pro	Pro	Lys	Pro	Lys	Asp
					245				250						255	
	Thr	Leu	Met	Ile	Ser	Arg	Thr	Pro	Glu	Val	Thr	Cys	Val	Val	Val	Asp
			260						265						270	
35	Val	Ser	His	Glu	Asp	Pro	Glu	Val	Lys	Phe	Asn	Trp	Tyr	Val	Asp	Gly
			275					280					285			
	Val	Glu	Val	His	Asn	Ala	Lys	Thr	Lys	Pro	Arg	Glu	Glu	Gln	Tyr	Asn
	290						295					300				
	Ser	Thr	Tyr	Arg	Val	Val	Ser	Val	Leu	Thr	Val	Leu	His	Gln	Asp	Trp
40	305					310					315					320
	Leu	Asn	Gly	Lys	Glu	Tyr	Lys	Cys	Lys	Val	Ser	Asn	Lys	Ala	Leu	Pro
					325					330					335	
	Ala	Pro	Ile	Glu	Lys	Thr	Ile	Ser	Lys	Ala	Lys	Gly	Gln	Pro	Arg	Glu
				340					345					350		
	Pro	Gln	Val	Tyr	Thr	Leu	Pro	Pro	Ser	Arg	Glu	Glu	Met	Thr	Lys	Asn
45			355						360					365		
	Gln	Val	Ser	Leu	Thr	Cys	Leu	Val	Lys	Gly	Phe	Tyr	Pro	Ser	Asp	Ile
	370						375					380				
	Ala	Val	Glu	Trp	Glu	Ser	Asn	Gly	Gln	Pro	Glu	Asn	Asn	Tyr	Lys	Thr
	385					390					395					400
50	Thr	Pro	Pro	Val	Leu	Asp	Ser	Asp	Gly	Ser	Phe	Phe	Leu	Tyr	Ser	Lys
					405					410					415	
	Leu	Thr	Val	Asp	Lys	Ser	Arg	Trp	Gln	Gln	Gly	Asn	Val	Phe	Ser	Cys
			420						425					430		
	Ser	Val	Met	His	Glu	Ala	Leu	His	Asn	His	Tyr	Thr	Gln	Lys	Ser	Leu
			435						440					445		
55	Ser	Leu	Ser	Pro	Gly	Lys										
	450															

<211> 454  
 <212> PRT  
 <213> Artificial Sequence

5 <220>  
 <223> Hu9D5 H5 heavy chain

<400> 34

10           Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
               1                           5                           10                           15  
               Ser Leu Arg Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
                           20                           25                           30  
 15           Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val  
                           35                           40                           45  
               Ala Glu Ile Ser Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
                           50                           55                           60  
               Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
                           65                           70                           75                           80  
 20           Leu Gln Met Asn Leu Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys  
                           85                           90                           95  
               Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
                           100                           105                           110  
               Val Trp Gly Gln Gly Thr Thr Val Thr Val Ser Ser Ala Ser Thr Lys  
                           115                           120                           125  
 25           Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys Ser Thr Ser Gly

30

35

40

45

50

55

		130				135				140						
	Gly	Thr	Ala	Ala	Leu	Gly	Cys	Leu	Val	Lys	Asp	Tyr	Phe	Pro	Glu	Pro
	145					150					155					160
5	Val	Thr	Val	Ser	Trp	Asn	Ser	Gly	Ala	Leu	Thr	Ser	Gly	Val	His	Thr
					165						170					175
	Phe	Pro	Ala	Val	Leu	Gln	Ser	Ser	Gly	Leu	Tyr	Ser	Leu	Ser	Ser	Val
					180						185					190
	Val	Thr	Val	Pro	Ser	Ser	Ser	Leu	Gly	Thr	Gln	Thr	Tyr	Ile	Cys	Asn
					195						200					205
10	Val	Asn	His	Lys	Pro	Ser	Asn	Thr	Lys	Val	Asp	Lys	Arg	Val	Glu	Pro
					210						215					220
	Lys	Ser	Cys	Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala	Pro	Glu
	225					230						235				240
	Leu	Leu	Gly	Gly	Pro	Ser	Val	Phe	Leu	Phe	Pro	Pro	Lys	Pro	Lys	Asp
15					245						250					255
	Thr	Leu	Met	Ile	Ser	Arg	Thr	Pro	Glu	Val	Thr	Cys	Val	Val	Val	Asp
					260						265					270
	Val	Ser	His	Glu	Asp	Pro	Glu	Val	Lys	Phe	Asn	Trp	Tyr	Val	Asp	Gly
					275						280					285
20	Val	Glu	Val	His	Asn	Ala	Lys	Thr	Lys	Pro	Arg	Glu	Glu	Gln	Tyr	Asn
					290						295					300
	Ser	Thr	Tyr	Arg	Val	Val	Ser	Val	Leu	Thr	Val	Leu	His	Gln	Asp	Trp
	305					310						315				320
	Leu	Asn	Gly	Lys	Glu	Tyr	Lys	Cys	Lys	Val	Ser	Asn	Lys	Ala	Leu	Pro
					325							330				335
25	Ala	Pro	Ile	Glu	Lys	Thr	Ile	Ser	Lys	Ala	Lys	Gly	Gln	Pro	Arg	Glu
					340						345					350
	Pro	Gln	Val	Tyr	Thr	Leu	Pro	Pro	Ser	Arg	Glu	Glu	Met	Thr	Lys	Asn
					355						360					365
	Gln	Val	Ser	Leu	Thr	Cys	Leu	Val	Lys	Gly	Phe	Tyr	Pro	Ser	Asp	Ile
30					370							375				380
	Ala	Val	Glu	Trp	Glu	Ser	Asn	Gly	Gln	Pro	Glu	Asn	Asn	Tyr	Lys	Thr
	385					390						395				400
	Thr	Pro	Pro	Val	Leu	Asp	Ser	Asp	Gly	Ser	Phe	Phe	Leu	Tyr	Ser	Lys
					405							410				415
35	Leu	Thr	Val	Asp	Lys	Ser	Arg	Trp	Gln	Gln	Gly	Asn	Val	Phe	Ser	Cys
					420							425				430
	Ser	Val	Met	His	Glu	Ala	Leu	His	Asn	His	Tyr	Thr	Gln	Lys	Ser	Leu
					435							440				445
	Ser	Leu	Ser	Pro	Gly	Lys										
					450											

40  
 <210> 35  
 <211> 219  
 <212> PRT  
 <213> Artificial Sequence

45  
 <220>  
 <223> Hu9D5 L1 light chain

50  
 <400> 35

55

```

Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
 1          5          10          15
Glu Pro Ala Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
 20          25          30
5  Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Lys Pro Gly Gln Ser
 35          40          45
Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
 50          55          60
10 Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
 65          70          75          80
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His

          85          90          95
15 Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
 100          105          110
Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu
 115          120          125
Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe
 130          135          140
20 Tyr Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln
 145          150          155          160
Ser Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser
 165          170          175
25 Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu
 180          185          190
Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser
 195          200          205
Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys
 210          215

```

```

30 <210> 36
    <211> 219
    <212> PRT
    <213> Artificial Sequence
35
    <220>
    <223> Hu9D5 L2 light chain
40
    <400> 36

```

45

50

55

```

    Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
    1          5          10          15
    Glu Pro Ala Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
    20          25          30
5   Asn Gly Asn Thr Tyr Leu Tyr Trp Tyr Leu Gln Lys Pro Gly Gln Ser
    35          40          45
    Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
    50          55          60
10  Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
    65          70          75          80
    Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
    85          90          95
    Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
    100         105         110
15  Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu
    115         120         125
    Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe
    130         135         140
    Tyr Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln
    145         150         155         160
20  Ser Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser
    165         170         175
    Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu
    180         185         190
25  Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser
    195         200         205
    Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys
    210         215

```

```

30  <210> 37
    <211> 219
    <212> PRT
    <213> Artificial Sequence

```

```

35  <220>
    <223> Hu9D5 L3 light chain

```

```

    <400> 37

```

40

45

50

55

```

Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
 1           5           10           15
Glu Pro Ala Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
 20           25           30
5  Asn Gly Asn Thr Tyr Leu Tyr Trp Tyr Leu Gln Lys Pro Gly Gln Ser
 35           40           45
Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
 50           55           60
Ser Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
10 65           70           75
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
 85           90           95
Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
15 100          105          110
Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu
115          120          125
Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe
130          135          140
Tyr Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln
20 145          150          155
Ser Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser
165          170          175
Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu
25 180          185          190
Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser
195          200          205
Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys
210          215

```

<210> 38

30 <211> 219

<212> PRT

<213> Artificial Sequence

<220>

35 <223> Hu9D5 L4 light chain

<400> 38

```

Asp Ile Val Met Thr Gln Ser Ala Pro Ser Leu Pro Val Thr Pro Gly
40  1           5           10           15
Glu Pro Val Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
 20           25           30
Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg Pro Gly Gln Ser
 35           40           45
45  Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
 50           55           60
Ser Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Arg Ile
65           70           75
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
 85           90           95
Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
50 100          105          110
Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu
115          120          125
Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe
55 130          135          140
Tyr Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln
145          150          155

```

Ser Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser  
 165 170 175  
 Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu  
 180 185 190  
 5 Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser  
 195 200 205  
 Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys  
 210 215

10 <210> 39  
 <211> 219  
 <212> PRT  
 <213> Artificial Sequence

15 <220>  
 <223> Hu9D5 L5 light chain

<400> 39

20 Asp Ile Val Met Thr Gln Ser Ala Pro Ser Leu Pro Val Thr Pro Gly  
 1 5 10 15  
 Glu Ser Val Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser  
 20 25 30  
 25 Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg Pro Gly Gln Ser  
 35 40 45  
 Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro  
 50 55 60  
 Ser Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Arg Ile  
 65 70 75 80  
 30 Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His  
 85 90 95  
 Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys  
 100 105 110  
 Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu  
 115 120 125  
 35 Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe  
 130 135 140  
 Tyr Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln  
 145 150 155 160  
 40 Ser Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser  
 165 170 175  
 Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu  
 180 185 190  
 Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser  
 195 200 205  
 45 Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys  
 210 215

<210> 40  
 <211> 429  
 50 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> m9D5VH with signal peptide

55 <400> 40

atggactttg ggctcagett gatttttcott gtccttqgtt taaaagggtgt cctgtgtgaa 60  
 gtgaagctgg tggagtctgg gggagggtta gtcagcctg gagggtcctt gaaactctcc 120  
 tgtgcagcct ctggattcac tttcagtagc tataccatgt ctggggttcg ccagactcca 180  
 gaaaaagggc tggagttggt cgcagaaatt agtaaatagt gtgataccac ctactatcca 240  
 5 gacactgtaa agggccgatt caccttctcc agagacaatg ccaagaacac cctgtacctg 300

caaatgagca gtctgaagtc tgaggacacg gccatgtatt actgtgcaag acattattac 360  
 tatgggtggt gctacggggg gtggttcttc gatgtctggg gcacagggac cacggtcacc 420  
 10 gtctctctcg 429

- <210> 41
- <211> 143
- <212> PRT
- 15 <213> Artificial Sequence

- <220>
- <223> m9D5VH with signal peptide

20 <400> 41

	Met	Asp	Phe	Gly	Leu	Ser	Leu	Ile	Phe	Leu	Val	Leu	Val	Leu	Lys	Gly
	1				5					10					15	
25	Val	Leu	Cys	Glu	Val	Lys	Leu	Val	Glu	Ser	Gly	Gly	Gly	Leu	Val	Gln
				20					25					30		
	Pro	Gly	Gly	Ser	Leu	Lys	Leu	Ser	Cys	Ala	Ala	Ser	Gly	Phe	Thr	Phe
				35				40					45			
	Ser	Ser	Tyr	Thr	Met	Ser	Trp	Val	Arg	Gln	Thr	Pro	Glu	Lys	Arg	Leu
30				50			55					60				
	Glu	Leu	Val	Ala	Glu	Ile	Ser	Asn	Ser	Gly	Asp	Thr	Thr	Tyr	Tyr	Pro
	65				70						75					80
	Asp	Thr	Val	Lys	Gly	Arg	Phe	Thr	Phe	Ser	Arg	Asp	Asn	Ala	Lys	Asn
				85						90					95	
	Thr	Leu	Tyr	Leu	Gln	Met	Ser	Ser	Leu	Lys	Ser	Glu	Asp	Thr	Ala	Met
35				100					105					110		
	Tyr	Tyr	Cys	Ala	Arg	His	Tyr	Tyr	Tyr	Gly	Gly	Gly	Tyr	Gly	Gly	Trp
				115					120					125		
	Phe	Phe	Asp	Val	Trp	Gly	Thr	Gly	Thr	Thr	Val	Thr	Val	Ser	Ser	
40				130			135					140				

- <210> 42
- <211> 396
- <212> DNA
- 45 <213> Artificial Sequence

- <220>
- <223> m9D5VL with signal peptide

50 <400> 42

atgaggtgcc tagctgagtt cctggggctg cttgtgctct ggatccctgg agccattggg 60  
 gatattgtga tgactcaggc tgcaccctct gtacctgtca ctctggaga gtcagtatcc 120  
 atctcctgca ggtctagtaa gactctcctg catagtaatg gcaacactta cttgtattgg 180  
 55 ttcttcgaga ggccaggcca gtctctctca ctctgatat atcgggtgtc caaccttggc 240  
 tcaggagtcc cagacaggtt cagtggcagt gggtcaggaa ctgctttcac actgagaatc 300  
 agtagagtgg aggctgagga tgtgggtgtt tattactgta tgcaacattt agaatatccg 360  
 ctcaagttcg gtgctgggac caagctggag ctgaaa 396

<210> 43  
 <211> 132  
 <212> PRT  
 <213> Artificial Sequence

5

<220>  
 <223> m9D5VL with signal peptide

<400> 43

10

```
Met Arg Cys Leu Ala Glu Phe Leu Gly Leu Leu Val Leu Trp Ile Pro
 1           5           10           15
Gly Ala Ile Gly Asp Ile Val Met Thr Gln Ala Ala Pro Ser Val Pro
                20                25                30
```

15

```
Val Thr Pro Gly Glu Ser Val Ser Ile Ser Cys Arg Ser Ser Lys Ser
                35                40                45
```

20

```
Leu Leu His Ser Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg
 50           55           60
```

```
Pro Gly Gln Ser Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala
 65           70           75           80
```

```
Ser Gly Val Pro Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe
 85           90           95
```

25

```
Thr Leu Arg Ile Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr
 100          105          110
```

```
Cys Met Gln His Leu Glu Tyr Pro Leu Thr Phe Gly Ala Gly Thr Lys
 115          120          125
```

```
Leu Glu Leu Lys
 130
```

30

<210> 44  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

35

<220>  
 <223> Hu9D5VHv1

<400> 44

40

```
gaggtgcagc tggtaggagtc cggcggcggc ctggtgcagc cggcggctc cctggcctg 60
tcttgccgcg cctccggett caccttctcc tctacacca tgtcctgggt ggcctaggcc 120
ccgggcaagg gctcggagct ggtggccgag atctccaact cggcggacac cactactac 180
ccgacacccg tgaagggcgg ctccaccttc tcccgagaca acgccaagaa ctccctgtac 240
ctgcagtcca actccctgaa ggccgaggac accgcctgt actactggc cggccactac 300
tactacggcg gcgctacgg cggctggttc ttgcagctgt gggccaggg caccacctg 360
accgtgtcct ca                                     372
```

45

<210> 45  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

50

<220>  
 <223> Hu9D5VHv2

<400> 45

55

```

gaggtgcagc  tggtaggagtc  cggcggggcgc  ctggtgcagc  cggcgggctc  cctgogcctg  60
tctgtgcgcg  cctccggctt  caccttctcc  tctacacca  tgtcctgggt  gggccaggcc  120
cccggcaagg  gcctggagct  ggtggccgag  atctccaact  cggcggacac  cacctactac  180
cccgacaccg  tgaagggccg  cttcaccttc  tcccgogaca  acgccaagaa  ctccctgtac  240
ctgcagtcca  acctgctgcg  cggcggaggac  accgocgtgt  actactgggc  cggccactac  300
tactacggcg  gcggtacgg  cggctggcttc  ttogacgtgt  ggggccaggg  caccaccgtg  360
accgtgtcct  ca  372

```

5  
 <210> 46  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

15  
 <220>  
 <223> Hu9D5VHv2b

<400> 46

```

gaggtgcaac  tggtaggagtc  tggcggggcgc  ttggtgcaac  ctggcggctc  cctgaagctg  60
tctgtgcccg  cctccggctt  caccttccgc  agctatacca  tgtcttgggt  gggccaaacc  120
cccgagaaga  ggtggagtt  ggtggctgag  attagtaata  gggcogatac  cacctactat  180
cccgataccg  tgaagggccg  cttcaccttc  tccagagata  atgctaagaa  taccctgtat  240

```

```

ctgcaaatga  gtagcctgaa  gtctgaggat  accgctatgt  attattgtgc  tagacattat  300
tattatggcg  gcggtatgg  cggctggcttc  ttogatgtgt  ggggccaggg  caccctggtc  360
accgtgtcct  ca  372

```

30  
 <210> 47  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

35  
 <220>  
 <223> Hu9D5VHv3

<400> 47

```

gaggtgcagc  tggtaggagtc  cggcggggcgc  ctggtgcagc  cggcgggctc  cctgogcctg  60
tctgtgcgcg  cctccggctt  caccttctcc  tctacacca  tgtcctgggt  gggccaggcc  120
cccggcaagg  gcctggagtg  ggtgtccgag  atctccaact  cggcggacac  cacctactac  180
cccgacaccg  tgaagggccg  cttcaccttc  tcccgogaca  acgccaagaa  ctccctgtac  240
ctgcagtcca  acctgctgcg  cggcggaggac  accgocgtgt  actactgggc  cggccactac  300
tactacggcg  gcggtacgg  cggctggcttc  ttogacgtgt  ggggccaggg  caccaccgtg  360
accgtgtcct  ca  372

```

45  
 <210> 48  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

50  
 <220>  
 <223> Hu9D5VHv3b

<400> 48

55

```

gaggtgcaac  ttggtggagtc  tggcgggoggc  ttggtgcaac  ctggcgggtc  cctgaagctg  60
tcctgtgccc  cctccggcctt  caccttcagc  agctatacca  tgtcttgggt  ggcgcaaac  120
cccgggaaga  ggctggagtt  ggtggctgag  attagtaata  ggggogatac  cacctactat  180
cccgataccg  tgaagggccc  cttcaccatt  tccagagata  atgctaagaa  tacctctgat  240
ctgcaaatga  gtagcctgaa  gtctgaggat  accgctatgt  attattgtgc  tagacattat  300
tattatggcg  gcgctatgg  cggctggttc  ttogatgtgt  ggggccaagg  cacctctgtc  360
accgtgtcct  ca  372

```

5  
 <210> 49  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

15  
 <220>  
 <223> Hu9D5VHv4

<400> 49

```

gaggtgcagc  ttggtggagtc  cggcggcggc  ctggtgcagc  cggcgggtc  cctgcgcctg  60
tcctgcccgc  cctccggcctt  caccttcctc  tctacaccca  tgtcctgggt  ggcgcaaggc  120
cccggcaagg  gcttggagtg  ggtgtccgag  atctccaaat  cggcgagcac  cacctactac  180
cccgcacccg  tgaagggccc  cttcaccate  tcccgagaca  acgccaagaa  ctccctgtac  240
ctgcagtcac  acctgctgcg  cggcggaggc  accgctgtgt  actactgcgc  ccgccaactc  300
tactacggcg  gcgctacgg  cggctggttc  ttogatgtgt  ggggccaagg  caccacccgtg  360
accgtgtcct  ca  372

```

20  
 <210> 50  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

30  
 <220>  
 <223> Hu9D5VHv4b

35  
 <400> 50

```

gaggtgcaac  ttggtggagtc  tggcgggoggc  ttggtgcaac  ctggcgggtc  cctgaagctg  60
tcctgtgccc  cctccggcctt  caccttcagc  agctatacca  tgtcttgggt  ggcgcaaac  120
cccgggaaga  ggctggagtt  ggtggctgag  attagtaata  ggggogatac  cacctactat  180
cccgataccg  tgaagggccc  cttcaccatt  tccagagata  atgctaagaa  tacctctgat  240
ctgcaaatga  gtagcctgaa  gtctgaggat  accgctatgt  attattgtgc  tagacattat  300
tattatggcg  gcgctatgg  cggctggttc  ttogatgtgt  ggggccaagg  cacctctgtc  360
accgtgtcct  ca  372

```

45  
 <210> 51  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

50  
 <220>  
 <223> Hu9D5VHv5

<400> 51

55

```

gagggtgcagc  tggtaggagtc  eggcggggggc  ctgggtgcagc  cgggagggtc  cctgogcctg  60
tcctgcgccc  cctccggctt  caccttctcc  tctacacca  tgtcctgggt  gggccagacc  120
cccgagaagc  gcctggagct  ggtggccgag  atctccaact  ccggcgacac  cactactac  180
cccgacaacc  tgaagggccg  cttcaccatc  tcccgogaca  acgccaagaa  caccctgtac  240
ctgcagatga  acctgctgcg  cggccgaggac  aocggcgtgt  actactgogc  ccgccactac  300
tactacggcg  gggctacgg  cggctgggtc  ttogacgtgt  ggggccaagg  caccaccgtg  360
accgtgtcct  ca

```

5  
10  
<210> 52  
<211> 336  
<212> DNA  
<213> Artificial Sequence

15  
<220>  
<223> Hu9D5VLv1

<400> 52

```

gacatcgtga  tgaccagtc  cccctgtcc  ctgcccgtga  ccccgggoga  gcccgctcc  60
atctcctgcc  gtcctccaa  gtccctgctg  cactccaacg  gcaacacctc  cctgtactgg  120
ttcctgcaga  agcccggcca  gtccccccag  ctgctgatct  accgctgttc  caacctggcc  180
tccggcgtgc  ccgaccgctt  ctccggctcc  ggtccgggca  ccgacttcac  cctgaagatc  240
tcccgcgtgg  aggcccaggc  cgtgggcgtg  tactactgca  tgcagcaact  ggagtacccc  300
ctgaccttcg  gccagggcac  caagctggag  atcaaa

```

20  
25  
<210> 53  
<211> 336  
<212> DNA  
<213> Artificial Sequence

30  
<220>  
<223> Hu9D5VLv2

<400> 53

```

gacatcgtga  tgaccagtc  cccctgtcc  ctgcccgtga  ccccgggoga  gcccgctcc  60
atctcctgcc  gtcctccaa  gtccctgctg  cactccaacg  gcaacacctc  cctgtactgg  120
tacctgcaga  agcccggcca  gtccccccag  ctgctgatct  accgctgttc  caacctggcc  180
tccggcgtgc  ccgaccgctt  ctccggctcc  ggtccgggca  ccgacttcac  cctgaagatc  240
tcccgcgtgg  aggcccaggc  cgtgggcgtg  tactactgca  tgcagcaact  ggagtacccc  300
ctgaccttcg  gccagggcac  caagctggag  atcaaa

```

35  
40  
45  
<210> 54  
<211> 336  
<212> DNA  
<213> Artificial Sequence

50  
<220>  
<223> Hu9D5VLv3

<400> 54

```

gacatcgtga  tgaccagtc  cccctgtcc  ctgcccgtga  ccccgggoga  gcccgctcc  60
atctcctgcc  gtcctccaa  gtccctgctg  cactccaacg  gcaacacctc  cctgtactgg  120
tacctgcaga  agcccggcca  gtccccccag  ctgctgatct  accgctgttc  caacctggcc  180
tccggcgtgc  cctcccgctt  ctccggctcc  ggtccgggca  ccgacttcac  cctgaagatc  240
tcccgcgtgg  aggcccaggc  cgtgggcgtg  tactactgca  tgcagcaact  ggagtacccc  300
ctgaccttcg  gccagggcac  caagctggag  atcaaa

```

5  
 <210> 55  
 <211> 336  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Hu9D5VL4

<400> 55

```

gacatcgtga tgacccagtc cgtccccctcc ctgcccgtga cccccggcga gcccggtgac 60
atctcctgcc gctcctccaa gtccctgctg cactccaaag gcaacacctc cctgtactgg 120
ttcctgcagc gccccggcca gtccccccag ctgctgatct acccggtgtc caacctggcc 180
tcggcgtgc cctcccgtt ctccggctcc ggctccggca cgccttcac cctgggcata 240
tcccgcgtgg aggcggagga cgtgggcgtg tactactgca tgcagcacct ggagtacccc 300
ctgaccttcg gccaaaggcac caagctggag atcaaaa 336
    
```

20  
 <210> 56  
 <211> 336  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Hu9D5VL5

<400> 56

```

gacatcgtga tgacccagtc cgtccccctcc ctgcccgtga cccccggcga gtcogtgtcc 60
atctcctgcc gctcctccaa gtccctgctg cactccaaag gcaacacctc cctgtactgg 120
ttcctgcagc gccccggcca gtccccccag ctgctgatct acccggtgtc caacctggcc 180
tcggcgtgc cctcccgtt ctccggctcc ggctccggca cgccttcac cctgggcata 240
tcccgcgtgg aggcggagga cgtgggcgtg tactactgca tgcagcacct ggagtacccc 300
ctgaccttcg gccaaaggcac caagctggag atcaaaa 336
    
```

35  
 <210> 57  
 <211> 19  
 <212> PRT  
 <213> Artificial Sequence

<220>  
 <223> 9D5VH Signal Peptide

<400> 57

```

Met Asp Phe Gly Leu Ser Leu Ile Phe Leu Val Leu Val Leu Lys Gly
  1           5           10           15
Val Leu Cys
    
```

50  
 <210> 58  
 <211> 57  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> 9D5VH Signal Peptide

<400> 58  
atggacttg ggctcagctt gattttcctt gtcttgttt taaaagggtg cctgtgt 57

5  
<210> 59  
<211> 20  
<212> PRT  
<213> Artificial Sequence

10  
<220>  
<223> 9D5VL Signal Peptide

<400> 59

15	Met	Arg	Cys	Leu	Ala	Glu	Phe	Leu	Gly	Leu	Leu	Val	Leu	Trp	Ile	Pro
	1				5				10					15		
	Gly	Ala	Ile	Gly												
				20												

20  
<210> 60  
<211> 60  
<212> DNA  
<213> Artificial Sequence

25  
<220>  
<223> 9D5VL Signal Peptide

<400> 60  
atgagggtgcc tagctgagtt cctggggctg cttgtgctct ggatccctgg agccattggg 60

30  
<210> 61  
<211> 124  
<212> PRT  
<213> Artificial Sequence

35  
<220>  
<223> m14G8VH

<400> 61

40	Glu	Val	Lys	Leu	Val	Glu	Ser	Gly	Gly	Gly	Leu	Val	Gln	Pro	Gly	Gly
	1				5				10						15	
	Ser	Leu	Lys	Leu	Ser	Cys	Ala	Ala	Ser	Gly	Phe	Thr	Phe	Ser	Ser	Tyr
				20					25					30		
	Thr	Met	Ser	Trp	Val	Arg	Gln	Thr	Pro	Glu	Lys	Arg	Leu	Glu	Leu	Val
45			35				40						45			
	Ala	Gln	Ile	Asn	Asn	Ser	Gly	Asp	Thr	Thr	Tyr	Tyr	Pro	Asp	Thr	Val
	50					55					60					
	Lys	Gly	Arg	Phe	Thr	Ile	Ser	Arg	Asp	Asn	Ala	Lys	Asn	Thr	Leu	Tyr
	65				70					75					80	
50	Leu	Gln	Met	Ser	Ser	Leu	Lys	Ser	Glu	Asp	Thr	Ala	Met	Tyr	Tyr	Cys
				85					90					95		
	Ala	Arg	His	Tyr	Tyr	Tyr	Gly	Gly	Gly	Tyr	Gly	Gly	Trp	Phe	Phe	Asp
			100					105						110		
	Val	Trp	Gly	Thr	Gly	Thr	Thr	Val	Thr	Val	Ser	Ser				
55			115					120								

<210> 62  
<211> 118  
<212> PRT

<213> Artificial Sequence

<220>

<223> 1MQK\_H

5

<400> 62

```

    Glu Val Lys Leu Gln Glu Ser Gly Gly Asp Leu Val Gln Pro Gly Gly
    1           5           10           15
    Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr
    20           25           30
    Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Trp Val
    35           40           45
    Ala Ser Ile Asn Asn Gly Gly Gly Arg Thr Tyr Tyr Pro Asp Thr Val
    50           55           60
    Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr
    65           70           75           80
    Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys
    85           90           95
    Val Arg His Gln Tyr Tyr Tyr Ala Met Asp Tyr Trp Gly Gln Gly Thr
    100          105          110
    Thr Val Thr Val Ser Ser
    115
  
```

<210> 63

25

<211> 124

<212> PRT

<213> Artificial Sequence

<220>

30

<223> AAD30410

<400> 63

```

    Gln Val Gln Leu Val Gln Ser Gly Gly Gly Leu Val Gln Pro Gly Gly
    1           5           10           15
    Ser Arg Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr
    20           25           30
    Ala Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Trp Val
    35           40           45
    Ala Ala Ile Ser Thr Asp Gly Ser Phe Ile Tyr Tyr Ala Asp Thr Val
    50           55           60
    Lys Gly Arg Phe Thr Ile Ser Gly Asp Asn Ser Lys Asn Thr Leu Tyr
    65           70           75           80
    Leu Gln Met Asn Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys
    85           90           95
    Ala Lys Asp Arg Gly Ile Asp Ala Thr Ala Gln Val Gly Arg Phe Asp
    100          105          110
    Pro Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser
    115          120
  
```

50

<210> 64

<211> 124

<212> PRT

<213> Artificial Sequence

55

<220>

<223> Hu14G8VHv1

<400> 64

Glu Val Lys Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
 1 5 10 15  
 Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
 5  
 Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val  
 20 25 30  
 Ala Glu Ile Asn Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 35 40 45  
 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 50 55 60  
 Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
 65 70 75 80  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 85 90 95  
 Val Trp Gly Thr Gly Thr Leu Val Thr Val Ser Ser  
 100 105 110  
 115 120

20 <210> 65  
 <211> 124  
 <212> PRT  
 <213> Artificial Sequence

25 <220>  
 <223> Hu14G8VHv2

<400> 65

Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
 1 5 10 15  
 Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
 20 25 30  
 Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val  
 35 40 45  
 Ala Glu Ile Asn Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 65 70 75 80  
 Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
 85 90 95  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 100 105 110  
 Val Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser  
 115 120

45 <210> 66  
 <211> 124  
 <212> PRT  
 <213> Artificial Sequence

50 <220>  
 <223> Hu14G8VHv3

<400> 66

55

```

    Glu Val Gln Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly
    1           5           10           15
    Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr
    20           25           30
5   Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val
    35           40           45
    Ala Glu Ile Asn Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val
    50           55           60
10  Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr
    65           70           75           80
    Leu Gln Met Asn Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys
    85           90           95

    Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp
15  Val Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser
    100           105           110
    115           120

```

20 <210> 67  
 <211> 5  
 <212> PRT  
 <213> Artificial Sequence

25 <220>  
 <223> 14G8 CDRH1 - Kabat  
 <400> 67

```

30          Ser Tyr Thr Met Ser
              1           5

```

35 <210> 68  
 <211> 17  
 <212> PRT  
 <213> Artificial Sequence

40 <220>  
 <223> 14G8 CDRH2 - Kabat  
 <400> 68

```

45  Glu Ile Asn Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val Lys
    1           5           10           15
    Gly

```

50 <210> 69  
 <211> 15  
 <212> PRT  
 <213> Artificial Sequence

55 <220>  
 <223> 14G8 CDRH3 - Kabat  
 <400> 69

```

    His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp Val
    1           5           10           15

```

<210> 70  
 <211> 113  
 <212> PRT  
 <213> Artificial Sequence

5

<220>  
 <223> m14G8VL

<400> 70

10

```

Asp Ile Val Met Thr Gln Ala Ala Pro Ser Val Pro Val Thr Pro Gly
 1           5           10           15
Glu Ser Val Ser Ile Ser Cys Arg Ser Asn Lys Ser Leu Leu His Ser
 20
Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg Pro Gly Gln Ser
 35           40           45
Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
    
```

15

20

```

          50           55           60
Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Arg Ile
 65           70           75           80
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
          85           90           95
Leu Glu Tyr Pro Leu Thr Phe Gly Ala Gly Thr Lys Leu Glu Leu Lys
 25           100          105          110
Arg
    
```

25

<210> 71  
 <211> 113  
 <212> PRT  
 <213> Artificial Sequence

30

<220>  
 <223> MJU\_L

35

<400> 71

```

Asp Ile Val Met Thr Gln Ala Ala Pro Ser Val Pro Val Thr Pro Gly
 1           5           10           15
Glu Ser Val Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
 20
Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg Pro Gly Gln Ser
 35           40           45
Pro Gln Leu Leu Ile Tyr Arg Met Ser Asn Leu Ala Ser Gly Val Pro
 50           55           60
Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Arg Ile
 65           70           75           80
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Leu Gln His
          85           90           95
Leu Glu Tyr Pro Phe Thr Phe Gly Ala Gly Thr Lys Leu Glu Leu Lys
 50           100          105          110
Arg
    
```

40

45

50

<210> 72  
 <211> 113  
 <212> PRT  
 <213> Artificial Sequence

55

<220>  
<223> ABA71374

<400> 72

5

10

15

20

```

Asp Ile Val Met Thr Gln Thr Pro Leu Ser Leu Pro Val Thr Pro Gly
 1           5           10           15
Glu Ser Ala Ser Ile Ser Cys Arg Ser Ser Gln Ser Leu Arg His Tyr
                20           25           30
Ser Gly Tyr Thr Tyr Ile Asp Trp Tyr Leu Gln Lys Pro Gly Gln Ser
                35           40           45
Pro Gln Val Leu Ile Tyr Leu Gly Ser Asn Arg Ala Ser Gly Val Pro
 50           55           60
Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
65           70           75           80
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln Ala
                85           90           95
Leu Gln Thr Val Cys Ser Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
                100           105           110
Arg
    
```

<210> 73  
<211> 113  
<212> PRT  
<213> Artificial Sequence

25

<220>  
<223> ABC66952

<400> 73

30

35

40

45

```

Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
 1           5           10           15
Gln Pro Ala Ser Ile Ser Cys Arg Ser Asn Gln Ser Leu Leu Tyr Ser
                20           25           30
Asn Gly Tyr Asn Tyr Leu Asp Trp Tyr Leu Gln Lys Pro Gly Gln Ser
                35           40           45
Pro Gln Leu Leu Ile Tyr Ser Gly Ser Asn Arg Ala Ser Gly Val Pro
 50           55           60
Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
65           70           75           80
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln Ala
                85           90           95
Leu Gln Ser Pro Tyr Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
                100           105           110
Arg
    
```

<210> 74  
<211> 113  
<212> PRT  
<213> Artificial Sequence

50

<220>  
<223> Hu14G8VLv1

55

<400> 74

```

    Asp Ile Val Met Thr Gln Ser Ala Pro Ser Leu Pro Val Thr Pro Gly
    1           5           10           15
    Glu Ser Val Ser Ile Ser Cys Arg Ser Asn Lys Ser Leu Leu His Ser
    20           25           30
5   Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Lys Pro Gly Gln Ser
    35           40           45
    Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
    50           55           60
10  Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Lys Ile
    65           70           75           80
    Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
    85           90           95
    Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
    100          105          110
15  Arg

```

```

<210> 75
<211> 113
<212> PRT
20  <213> Artificial Sequence

<220>
<223> Hu14G8VLv2

25  <400> 75

```

```

    Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
    1           5           10           15
30  Glu Pro Ala Ser Ile Ser Cys Arg Ser Asn Lys Ser Leu Leu His Ser
    20           25           30
    Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Lys Pro Gly Gln Ser
    35           40           45
    Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
    50           55           60
35  Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
    65           70           75           80
    Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
    85           90           95
    Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
    100          105          110
40  Arg

```

```

<210> 76
<211> 113
45  <212> PRT
    <213> Artificial Sequence

<220>
<223> Hu14G8VLv3

50  <400> 76

```

55

```

    Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
    1           5           10           15
    Glu Pro Ala Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
    20           25           30
5   Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Lys Pro Gly Gln Ser
    35           40           45
    Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
    50           55           60
10  Ser Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
    65           70           75           80
    Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
    85           90           95
    Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
    100          105          110
15  Arg

```

```

<210> 77
<211> 16
<212> PRT
20 <213> Artificial Sequence

<220>
<223> 14G8 CDRL1 (mouse, HuVLv1, HuVLv2) - Kabat

```

```

25 <400> 77

    Arg Ser Asn Lys Ser Leu Leu His Ser Asn Gly Asn Thr Tyr Leu Tyr
    1           5           10           15

```

```

30 <210> 78
<211> 7
<212> PRT
<213> Artificial Sequence

```

```

35 <220>
<223> 14G8 CDRL2 - Kabat

<400> 78

```

```

40

    Arg Val Ser Asn Leu Ala Ser
    1           5

```

```

45 <210> 79
<211> 9
<212> PRT
<213> Artificial Sequence

```

```

50 <220>
<223> 14G8 CDRL3 - Kabat

<400> 79

```

```

55

    Met Gln His Leu Glu Tyr Pro Leu Thr
    1           5

```

```

<210> 80
<211> 16

```

<212> PRT  
 <213> Artificial Sequence

<220>  
 <223> 14G8 CDRL1 (HuVLv3) - Kabat

<400> 80

10 Arg Ser Ser Lys Ser Leu Leu His Ser Asn Gly Asn Thr Tyr Leu Tyr  
 1 5 10 15

<210> 81  
 <211> 454  
 <212> PRT  
 <213> Artificial Sequence

<220>  
 <223> Hu14G8 H1 heavy chain

<400> 81

25 Glu Val Lys Leu Val Glu Ser Gly Gly Gly Leu Val Gln Pro Gly Gly  
 1 5 10 15  
 Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe Ser Ser Tyr  
 20 25 30  
 Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu Glu Leu Val  
 35 40 45  
 Ala Glu Ile Asn Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 30 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 65 70 75 80  
 Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met Tyr Tyr Cys  
 85 90 95  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 100 105 110  
 35 Val Trp Gly Thr Gly Thr Leu Val Thr Val Ser Ser Ala Ser Thr Lys  
 115 120 125  
 Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys Ser Thr Ser Gly  
 130 135 140

40

45

50

55

Gly Thr Ala Ala Leu Gly Cys Leu Val Lys Asp Tyr Phe Pro Glu Pro  
 145 150 155 160  
 Val Thr Val Ser Trp Asn Ser Gly Ala Leu Thr Ser Gly Val His Thr  
 165 170 175  
 5 Phe Pro Ala Val Leu Gln Ser Ser Gly Leu Tyr Ser Leu Ser Ser Val  
 180 185 190  
 Val Thr Val Pro Ser Ser Ser Leu Gly Thr Gln Thr Tyr Ile Cys Asn  
 195 200 205  
 10 Val Asn His Lys Pro Ser Asn Thr Lys Val Asp Lys Arg Val Glu Pro  
 210 215 220  
 Lys Ser Cys Asp Lys Thr His Thr Cys Pro Pro Cys Pro Ala Pro Glu  
 225 230 235 240  
 Leu Leu Gly Gly Pro Ser Val Phe Leu Phe Pro Pro Lys Pro Lys Asp  
 245 250 255  
 15 Thr Leu Met Ile Ser Arg Thr Pro Glu Val Thr Cys Val Val Val Asp  
 260 265 270  
 Val Ser His Glu Asp Pro Glu Val Lys Phe Asn Trp Tyr Val Asp Gly  
 275 280 285  
 Val Glu Val His Asn Ala Lys Thr Lys Pro Arg Glu Glu Gln Tyr Asn  
 290 295 300  
 20 Ser Thr Tyr Arg Val Val Ser Val Leu Thr Val Leu His Gln Asp Trp  
 305 310 315 320  
 Leu Asn Gly Lys Glu Tyr Lys Cys Lys Val Ser Asn Lys Ala Leu Pro  
 325 330 335  
 25 Ala Pro Ile Glu Lys Thr Ile Ser Lys Ala Lys Gly Gln Pro Arg Glu  
 340 345 350  
 Pro Gln Val Tyr Thr Leu Pro Pro Ser Arg Glu Glu Met Thr Lys Asn  
 355 360 365  
 Gln Val Ser Leu Thr Cys Leu Val Lys Gly Phe Tyr Pro Ser Asp Ile  
 370 375 380  
 30 Ala Val Glu Trp Glu Ser Asn Gly Gln Pro Glu Asn Asn Tyr Lys Thr  
 385 390 395 400  
 Thr Pro Pro Val Leu Asp Ser Asp Gly Ser Phe Phe Leu Tyr Ser Lys  
 405 410 415  
 Leu Thr Val Asp Lys Ser Arg Trp Gln Gln Gly Asn Val Phe Ser Cys  
 420 425 430  
 35 Ser Val Met His Glu Ala Leu His Asn His Tyr Thr Gln Lys Ser Leu  
 435 440 445  
 Ser Leu Ser Pro Gly Lys  
 450

40 <210> 82  
 <211> 454  
 <212> PRT  
 <213> Artificial Sequence

45 <220>  
 <223> Hu14G8 H2 heavy chain

<400> 82

50

55

	Glu	Val	Gln	Leu	Val	Glu	Ser	Gly	Gly	Gly	Leu	Val	Gln	Pro	Gly	Gly
	1			5					10						15	
	Ser	Leu	Lys	Leu	Ser	Cys	Ala	Ala	Ser	Gly	Phe	Thr	Phe	Ser	Ser	Tyr
5				20					25					30		
	Thr	Met	Ser	Trp	Val	Arg	Gln	Thr	Pro	Glu	Lys	Arg	Leu	Glu	Leu	Val
				35				40					45			
	Ala	Glu	Ile	Asn	Asn	Ser	Gly	Asp	Thr	Thr	Tyr	Tyr	Pro	Asp	Thr	Val
	50						55					60				
10	Lys	Gly	Arg	Phe	Thr	Ile	Ser	Arg	Asp	Asn	Ala	Lys	Asn	Thr	Leu	Tyr
	65					70					75					80
	Leu	Gln	Met	Ser	Ser	Leu	Lys	Ser	Glu	Asp	Thr	Ala	Met	Tyr	Tyr	Cys
				85					90						95	
15	Ala	Arg	His	Tyr	Tyr	Tyr	Gly	Gly	Gly	Tyr	Gly	Gly	Trp	Phe	Phe	Asp
				100					105					110		
	Val	Trp	Gly	Gln	Gly	Thr	Leu	Val	Thr	Val	Ser	Ser	Ala	Ser	Thr	Lys
				115				120					125			
	Gly	Pro	Ser	Val	Phe	Pro	Leu	Ala	Pro	Ser	Ser	Lys	Ser	Thr	Ser	Gly
20							135					140				
	Gly	Thr	Ala	Ala	Leu	Gly	Cys	Leu	Val	Lys	Asp	Tyr	Phe	Pro	Glu	Pro
	145					150					155					160
	Val	Thr	Val	Ser	Trp	Asn	Ser	Gly	Ala	Leu	Thr	Ser	Gly	Val	His	Thr
				165					170						175	
	Phe	Pro	Ala	Val	Leu	Gln	Ser	Ser	Gly	Leu	Tyr	Ser	Leu	Ser	Ser	Val
25				180					185					190		
	Val	Thr	Val	Pro	Ser	Ser	Ser	Leu	Gly	Thr	Gln	Thr	Tyr	Ile	Cys	Asn
				195				200					205			
	Val	Asn	His	Lys	Pro	Ser	Asn	Thr	Lys	Val	Asp	Lys	Arg	Val	Glu	Pro
				210			215				220					
30	Lys	Ser	Cys	Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala	Pro	Glu
	225					230					235					240
	Leu	Leu	Gly	Gly	Pro	Ser	Val	Phe	Leu	Phe	Pro	Pro	Lys	Pro	Lys	Asp
				245					250						255	
	Thr	Leu	Met	Ile	Ser	Arg	Thr	Pro	Glu	Val	Thr	Cys	Val	Val	Val	Asp
				260					265						270	
35	Val	Ser	His	Glu	Asp	Pro	Glu	Val	Lys	Phe	Asn	Trp	Tyr	Val	Asp	Gly
				275				280					285			
	Val	Glu	Val	His	Asn	Ala	Lys	Thr	Lys	Pro	Arg	Glu	Glu	Gln	Tyr	Asn
				290			295					300				
40	Ser	Thr	Tyr	Arg	Val	Val	Ser	Val	Leu	Thr	Val	Leu	His	Gln	Asp	Trp
	305					310					315					320
	Leu	Asn	Gly	Lys	Glu	Tyr	Lys	Cys	Lys	Val	Ser	Asn	Lys	Ala	Leu	Pro
				325						330					335	
	Ala	Pro	Ile	Glu	Lys	Thr	Ile	Ser	Lys	Ala	Lys	Gly	Gln	Pro	Arg	Glu
				340					345					350		
45	Pro	Gln	Val	Tyr	Thr	Leu	Pro	Pro	Ser	Arg	Glu	Glu	Met	Thr	Lys	Asn
				355				360					365			
	Gln	Val	Ser	Leu	Thr	Cys	Leu	Val	Lys	Gly	Phe	Tyr	Pro	Ser	Asp	Ile
				370			375					380				
	Ala	Val	Gln	Trp	Glu	Ser	Asn	Gly	Gln	Pro	Glu	Asn	Asn	Tyr	Lys	Thr
				385			390				395					400
50	Thr	Pro	Pro	Val	Leu	Asp	Ser	Asp	Gly	Ser	Phe	Phe	Leu	Tyr	Ser	Lys
				405						410					415	
	Leu	Thr	Val	Asp	Lys	Ser	Arg	Trp	Gln	Gln	Gly	Asn	Val	Phe	Ser	Cys
				420					425					430		
	Ser	Val	Met	His	Glu	Ala	Leu	His	Asn	His	Tyr	Thr	Gln	Lys	Ser	Leu
				435				440					445			
55	Ser	Leu	Ser	Pro	Gly	Lys										
				450												

<211> 454  
 <212> PRT  
 <213> Artificial Sequence

5 <220>  
 <223> Hu14G8 H3 heavy chain

<400> 83

10	Gln	Val	Gln	Leu	Val	Glu	Ser	Gly	Gly	Gly	Leu	Val	Gln	Pro	Gly	Gly
	1				5					10					15	
	Ser	Leu	Lys	Leu	Ser	Cys	Ala	Ala	Ser	Gly	Phe	Thr	Phe	Ser	Ser	Tyr
			20						25					30		
15	Thr	Met	Ser	Trp	Val	Arg	Gln	Thr	Pro	Glu	Lys	Arg	Leu	Glu	Leu	Val
			35					40					45			

20

25

30

35

40

45

50

55

Ala Glu Ile Asn Asn Ser Gly Asp Thr Thr Tyr Tyr Pro Asp Thr Val  
 50 55 60  
 Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn Thr Leu Tyr  
 65 70 75 80  
 5 Leu Gln Met Asn Ser Leu Arg Ala Glu Asp Thr Ala Val Tyr Tyr Cys  
 85 90 95  
 Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp Phe Phe Asp  
 100 105 110  
 10 Val Trp Gly Gln Gly Thr Leu Val Thr Val Ser Ser Ala Ser Thr Lys  
 115 120 125  
 Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys Ser Thr Ser Gly  
 130 135 140  
 Gly Thr Ala Ala Leu Gly Cys Leu Val Lys Asp Tyr Phe Pro Glu Pro  
 145 150 155 160  
 15 Val Thr Val Ser Trp Asn Ser Gly Ala Leu Thr Ser Gly Val His Thr  
 165 170 175  
 Phe Pro Ala Val Leu Gln Ser Ser Gly Leu Tyr Ser Leu Ser Ser Val  
 180 185 190  
 Val Thr Val Pro Ser Ser Ser Leu Gly Thr Gln Thr Tyr Ile Cys Asn  
 195 200 205  
 20 Val Asn His Lys Pro Ser Asn Thr Lys Val Asp Lys Arg Val Glu Pro  
 210 215 220  
 Lys Ser Cys Asp Lys Thr His Thr Cys Pro Pro Cys Pro Ala Pro Glu  
 225 230 235 240  
 25 Leu Leu Gly Gly Pro Ser Val Phe Leu Phe Pro Pro Lys Pro Lys Asp  
 245 250 255  
 Thr Leu Met Ile Ser Arg Thr Pro Gln Val Thr Cys Val Val Val Asp  
 260 265 270  
 Val Ser His Glu Asp Pro Glu Val Lys Phe Asn Trp Tyr Val Asp Gly  
 275 280 285  
 30 Val Glu Val His Asn Ala Lys Thr Lys Pro Arg Glu Glu Gln Tyr Asn  
 290 295 300  
 Ser Thr Tyr Arg Val Val Ser Val Leu Thr Val Leu His Gln Asp Trp  
 305 310 315 320  
 Leu Asn Gly Lys Glu Tyr Lys Cys Lys Val Ser Asn Lys Ala Leu Pro  
 325 330 335  
 35 Ala Pro Ile Glu Lys Thr Ile Ser Lys Ala Lys Gly Gln Pro Arg Glu  
 340 345 350  
 Pro Gln Val Tyr Thr Leu Pro Pro Ser Arg Glu Glu Met Thr Lys Asn  
 355 360 365  
 40 Gln Val Ser Leu Thr Cys Leu Val Lys Gly Phe Tyr Pro Ser Asp Ile  
 370 375 380  
 Ala Val Glu Trp Glu Ser Asn Gly Gln Pro Glu Asn Asn Tyr Lys Thr  
 385 390 395 400  
 Thr Pro Pro Val Leu Asp Ser Asp Gly Ser Phe Phe Leu Tyr Ser Lys  
 405 410 415  
 45 Leu Thr Val Asp Lys Ser Arg Trp Gln Gln Gly Asn Val Phe Ser Cys  
 420 425 430  
 Ser Val Met His Glu Ala Leu His Asn His Tyr Thr Gln Lys Ser Leu  
 435 440 445  
 Ser Leu Ser Pro Gly Lys  
 450

50

<210> 84  
 <211> 219  
 <212> PRT  
 <213> Artificial Sequence

55

<220>  
 <223> Hu14G8 L1 light chain

<400> 84

```

5      Asp Ile Val Met Thr Gln Ser Ala Pro Ser Leu Pro Val Thr Pro Gly
      1      5      10      15
      Glu Ser Val Ser Ile Ser Cys Arg Ser Asn Lys Ser Leu Leu His Ser
      20      35      40      45
      Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Lys Pro Gly Gln Ser
      50      55      60
      Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
      65      70      75      80
      Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe Thr Leu Lys Ile
      85      90      95
      Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
      100      105      110
      Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
      115      120      125
      Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu
      130      135      140
      Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe
      145      150      155      160
      Tyr Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln
      165      170      175
      Ser Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser
      180      185      190
      Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu
      195      200      205
      Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser
      210      215
      Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys

```

30

<210> 85  
 <211> 219  
 <212> PRT  
 <213> Artificial Sequence

35

<220>  
 <223> Hu14G8 L2 light chain

40

<400> 85

45

50

55

Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly  
 1 5 10 15  
 Glu Pro Ala Ser Ile Ser Cys Arg Ser Asn Lys Ser Leu Leu His Ser  
 20 25 30  
 5 Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Lys Pro Gly Gln Ser  
 35 40 45  
 Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro  
 50 55 60  
 10 Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile  
 65 70 75 80  
 Ser Arg Val Gln Ala Gln Asp Val Gly Val Tyr Tyr Cys Met Gln His  
 85 90 95  
 Leu Gln Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Gln Ile Lys  
 100 105 110  
 15 Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu  
 115 120 125  
 Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe  
 130 135 140  
 Tyr Pro Arg Gln Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln  
 145 150 155 160  
 20 Ser Gly Asn Ser Gln Gln Ser Val Thr Gln Gln Asp Ser Lys Asp Ser  
 165 170 175  
 Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu  
 180 185 190  
 25 Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser  
 195 200 205  
 Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys  
 210 215

30 <210> 86  
 <211> 219  
 <212> PRT  
 <213> Artificial Sequence

35 <220>  
 <223> Hu14G8 L3 light chain

<400> 86

40

45

50

55

```

Asp Ile Val Met Thr Gln Ser Pro Leu Ser Leu Pro Val Thr Pro Gly
 1           5           10           15
Glu Pro Ala Ser Ile Ser Cys Arg Ser Ser Lys Ser Leu Leu His Ser
 20           25           30
5  Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Lys Pro Gly Gln Ser
 35           40           45
Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala Ser Gly Val Pro
 50           55           60
Ser Arg Phe Ser Gly Ser Gly Ser Gly Thr Asp Phe Thr Leu Lys Ile
10 65           70           75
Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr Cys Met Gln His
 85           90           95
Leu Glu Tyr Pro Leu Thr Phe Gly Gln Gly Thr Lys Leu Glu Ile Lys
15 100          105          110
Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu
115          120          125
Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe
130          135          140
Tyr Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln
20 145          150          155
Ser Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser
165          170          175
Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu
180          185          190
Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser
25 195          200          205
Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys
 210          215

```

```

30 <210> 87
    <211> 429
    <212> DNA
    <213> Artificial Sequence

```

```

35 <220>
    <223> m14G8VH with signal peptide

```

<400> 87

```

40 atgaatttgc ggcacagctt gattttcctt gtccttggtt taaaagggtg cctgtgtgaa 60
gtgaagctgg tggagtctgg gggaggttta gtgcagcctg gagggtcctt gaaactctcc 120
tgtgcagcct ccggattcac ttccagttagc tataccatgt ottgggttgc ccagactcca 180
gaaaagagggc tggagttggt ccagaaaatt aataatagtg gtgataccac ctactatcca 240
gacactgtaa agggccgatt caccctctcc agagacaatg ccaagaacac cctgtacctg 300
caaatgagca gtctgaagtc tgaggacacg gccatgtatt actgtgcaag acattattac 360
45 tatggtggtg gctacggggg gtggttcttc gatgtctggg gcacagggac cacggtcacc 420
gtctctca 429

```

```

50 <210> 88
    <211> 143
    <212> PRT
    <213> Artificial Sequence

```

```

55 <220>
    <223> m14G8VH with signal peptide

```

<400> 88

```

Met Asn Phe Gly Leu Ser Leu Ile Phe Leu Val Leu Val Leu Lys Gly
 1          5          10          15
Val Leu Cys Glu Val Lys Leu Val Glu Ser Gly Gly Gly Leu Val Gln
 20          25          30
5 Pro Gly Gly Ser Leu Lys Leu Ser Cys Ala Ala Ser Gly Phe Thr Phe
 35          40          45
Ser Ser Tyr Thr Met Ser Trp Val Arg Gln Thr Pro Glu Lys Arg Leu
 50          55          60
Glu Leu Val Ala Glu Ile Asn Asn Ser Gly Asp Thr Thr Tyr Tyr Pro
 65          70          75          80
10 Asp Thr Val Lys Gly Arg Phe Thr Ile Ser Arg Asp Asn Ala Lys Asn
 85          90          95
Thr Leu Tyr Leu Gln Met Ser Ser Leu Lys Ser Glu Asp Thr Ala Met
 100          105          110
15 Tyr Tyr Cys Ala Arg His Tyr Tyr Tyr Gly Gly Gly Tyr Gly Gly Trp
 115          120          125
Phe Phe Asp Val Trp Gly Thr Gly Thr Thr Val Thr Val Ser Ser
 130          135          140

```

```

<210> 89
20 <211> 399
    <212> DNA
    <213> Artificial Sequence

<220>
25 <223> m14G8VL with signal peptide

<400> 89

```

```

30 atgagggtgcc tagctgagtt cctggggctg cttgtgctct ggatccctgg agccattggg 60
   gatattgtga tgactcaggc tgcaccctct gtacctgtca ctctgggaga gtcagtatcc 120
   atctcctgca ggtctaataa gagtctctct catagtaatg gcaacactta cttgtattgg 180
   ttctctgcaga ggcaggcca gtctctctcaa ctctgatat atcgggtgtc caaccttgcc 240
   tcaggagttc cagacaggtt cagtggcagt gggtcaggaa ctgctttcac actgagaatc 300
   agtagagttg aggctgagga tgtgggtgtt tattactgta tgcaacattt agaatatccg 360
35 ctacagttcg gtgctgggac caagctggag ctgaaacct 399

```

```

<210> 90
<211> 133
<212> PRT
40 <213> Artificial Sequence

<220>
<223> m14G8VL with signal peptide

45 <400> 90

```

```

Met Arg Cys Leu Ala Glu Phe Leu Gly Leu Leu Val Leu Trp Ile Pro
 1          5          10          15
Gly Ala Ile Gly Asp Ile Val Met Thr Gln Ala Ala Pro Ser Val Pro
 20          25          30
50 Val Thr Pro Gly Glu Ser Val Ser Ile Ser Cys Arg Ser Asn Lys Ser
 35          40          45
Leu Leu His Ser Asn Gly Asn Thr Tyr Leu Tyr Trp Phe Leu Gln Arg
 50          55          60
55

```

Pro Gly Gln Ser Pro Gln Leu Leu Ile Tyr Arg Val Ser Asn Leu Ala  
 65 70 75 80  
 Ser Gly Val Pro Asp Arg Phe Ser Gly Ser Gly Ser Gly Thr Ala Phe  
 85 90 95  
 5 Thr Leu Arg Ile Ser Arg Val Glu Ala Glu Asp Val Gly Val Tyr Tyr  
 100 105 110  
 Cys Met Gln His Leu Glu Tyr Pro Leu Thr Phe Gly Ala Gly Thr Lys  
 115 120 125  
 10 Leu Glu Leu Lys Arg  
 130

<210> 91  
 <211> 372  
 <212> DNA  
 15 <213> Artificial Sequence

<220>  
 <223> Hu14G8VHv1

20 <400> 91

gaggtgaagc tggtagagtc tggggggggc ttgggtgcaac ctggcgggctc cctgaagctg 60  
 tcctgtgccc cctccggcctt caccttcagc agctatacca tgtcttgggt ggcgcaaac 120  
 25 cccgagaaga ggctggagtt ggtggctgag attaataata gggcgatagc cacctactat 180  
 cccgataccg tgaagggccg cttcaccatt tccagagata atgctaagaa tacctctgat 240  
 ctgcaaatga gtagcctgaa gtctgaggat accgctatgt attattgtgc tagacattat 300  
 tattatggcg gggcctatgg cggctggctt ttogatgtgt ggggcaccgg cacctctgctc 360  
 accgtgtcct ca 372

30 <210> 92  
 <211> 372  
 <212> DNA  
 <213> Artificial Sequence

35 <220>  
 <223> Hu14G8VHv2

<400> 92

gaggtgcaac tggtagagtc tggggggggc ttgggtgcaac ctggcgggctc cctgaagctg 60  
 tcctgtgccc cctccggcctt caccttcagc agctatacca tgtcttgggt ggcgcaaac 120  
 40 cccgagaaga ggctggagtt ggtggctgag attaataata gggcgatagc cacctactat 180  
 cccgataccg tgaagggccg cttcaccatt tccagagata atgctaagaa tacctctgat 240  
 ctgcaaatga gtagcctgaa gtctgaggat accgctatgt attattgtgc tagacattat 300  
 45 tattatggcg gggcctatgg cggctggctt ttogatgtgt ggggcaccgg cacctctgctc 360  
 accgtgtcct ca 372

<210> 93  
 <211> 372  
 50 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> Hu14G8VHv3

55 <400> 93

```

gaggtgcaac  ttggtggagtc  tggcgggggc  ttggtgcaac  ctggcgggctc  cctgaagctg  60
tcctgtgccg  cctccgggctt  caccttcagc  agctatacca  tgtcttgggt  ggcgcaaac  120
cccgagaaga  ggctggagtt  ggtggctgag  attaataata  ggggogatac  cacctactat  180
cccgataccg  tgaaggggccg  cttcaccatt  tccagagata  atgctaagaa  taccctgtat  240
ctgcaaatga  atagcctgag  ggcctgaggat  accgctgtgt  attattgtgc  tagacattat  300
tattatggcg  gggctatgg  cggctggttc  ttcgatgtgt  ggggccaagg  cacctggctc  360
accgtgtcct  ca

```

5  
10  
<210> 94  
<211> 339  
<212> DNA  
<213> Artificial Sequence

15  
<220>  
<223> Hu14G8VLv1

<400> 94

```

gatatcgtga  tgaccocagtc  tggccccctcc  ctgcctgtga  cccctgggoga  gtcogtctcc  60
atctcctgcc  ggtcccaacaa  gagcctgctg  cacagcaacg  gcaacaacct  cctgtactgg  120
ttcctgcaaa  agcccgggcca  atccccctcaa  ctgctgatct  accgggtgtc  caacctggcc  180
tcgggctgtc  ccgataggtt  ctccggaaagc  ggctccggca  ccgcttccac  cctgaagatt  240
agttagagtcg  aggccgagga  tgtgggctgt  tactactgta  tgcaaacactt  ggagtaaccc  300
ctgacgttcg  gccaaaggcac  caagctggag  atcaagcgt

```

20  
25  
<210> 95  
<211> 339  
<212> DNA  
<213> Artificial Sequence

30  
<220>  
<223> Hu14G8VLv2

35  
<400> 95

```

gatatcgtga  tgaccocagtc  tccccctgtcc  ctgcctgtga  cccctgggoga  gcccgctccc  60
atctcctgcc  ggtcccaacaa  gagcctgctg  cacagcaacg  gcaacaacct  cctgtactgg  120
ttcctgcaaa  agcccgggcca  atccccctcaa  ctgctgatct  accgggtgtc  caacctggcc  180
tcgggctgtc  ccgataggtt  ctccggaaagc  ggctccggca  ccgatttccac  cctgaagatt  240
agttagagtcg  aggccgagga  tgtgggctgt  tactactgta  tgcaaacactt  ggagtaaccc  300
ctgacgttcg  gccaaaggcac  caagctggag  atcaagcgt

```

40  
45  
<210> 96  
<211> 339  
<212> DNA  
<213> Artificial Sequence

50  
<220>  
<223> Hu14G8VLv3

<400> 96

```

gatatcgtga  tgaccocagtc  tccccctgtcc  ctgcctgtga  cccctgggoga  gcccgctccc  60
atctcctgcc  ggtccagcaa  gagcctgctg  cacagcaacg  gcaacaacct  cctgtactgg  120
ttcctgcaaa  agcccgggcca  atccccctcaa  ctgctgatct  accgggtgtc  caacctggcc  180
tcgggctgtc  ccgataggtt  ctccggaaagc  ggctccggca  ccgatttccac  cctgaagatt  240
agttagagtcg  aggccgagga  tgtgggctgt  tactactgta  tgcaaacactt  ggagtaaccc  300
ctgacgttcg  gccaaaggcac  caagctggag  atcaagcgt

```

<210> 97  
 <211> 19  
 <212> PRT  
 <213> Artificial Sequence  
 5  
 <220>  
 <223> 14G8VH Signal Peptide  
 <400> 97  
 10  
     Met Asn Phe Gly Leu Ser Leu Ile Phe Leu Val Leu Val Leu Lys Gly  
         1                          5                          10                          15  
     Val Leu Cys  
 15  
 <210> 98  
 <211> 57  
 <212> DNA  
 <213> Artificial Sequence  
 20  
 <220>  
 <223> 14G8VH Signal Peptide  
 <400> 98  
 25  
 atgaattcg gctgagcgt gatttcctg gtgctggtgt tgaagggcgt gctgtgt 57  
 <210> 99  
 <211> 20  
 <212> PRT  
 <213> Artificial Sequence  
 30  
 <220>  
 <223> 14G8VL Signal Peptide  
 35  
 <400> 99  
     Met Arg Cys Leu Ala Glu Phe Leu Gly Leu Leu Val Leu Trp Ile Pro  
         1                          5                          10                          15  
     Gly Ala Ile Gly  
                           20  
 40  
 <210> 100  
 <211> 60  
 <212> DNA  
 <213> Artificial Sequence  
 45  
 <220>  
 <223> 14G8VL Signal Peptide  
 50  
 <400> 100  
 atgaggtagc tggccgagtt cctgggacct ctgggtcgtt ggatccctgg cgccatcggc 60  
 <210> 101  
 <211> 330  
 <212> PRT  
 <213> Artificial Sequence  
 55  
 <220>

&lt;223&gt; Human IgG1 heavy chain constant region

&lt;400&gt; 101

```

5      Ala Ser Thr Lys Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys
      1      5      10
      Ser Thr Ser Gly Gly Thr Ala Ala Leu Gly Cys Leu Val Lys Asp Tyr
      20      25      30
10     Phe Pro Glu Pro Val Thr Val Ser Trp Asn Ser Gly Ala Leu Thr Ser
      35      40      45
      Gly Val His Thr Phe Pro Ala Val Leu Gln Ser Ser Gly Leu Tyr Ser
      50      55      60
      Leu Ser Ser Val Val Thr Val Pro Ser Ser Ser Leu Gly Thr Gln Thr
      65      70      75      80
15     Tyr Ile Cys Asn Val Asn His Lys Pro Ser Asn Thr Lys Val Asp Lys
      85      90      95
      Arg Val Glu Pro Lys Ser Cys Asp Lys Thr His Thr Cys Pro Pro Cys
      100     105     110
      Pro Ala Pro Glu Leu Leu Gly Gly Pro Ser Val Phe Leu Phe Pro Pro
      115     120     125
20
      Lys Pro Lys Asp Thr Leu Met Ile Ser Arg Thr Pro Glu Val Thr Cys
      130     135     140
      Val Val Val Asp Val Ser His Glu Asp Pro Glu Val Lys Phe Asn Trp
      145     150     155     160
25     Tyr Val Asp Gly Val Glu Val His Asn Val Lys Thr Lys Pro Arg Glu
      165     170     175
      Glu Gln Tyr Asn Ser Thr Tyr Arg Val Val Ser Val Leu Thr Val Leu
      180     185     190
      His Gln Asp Trp Leu Asn Gly Lys Glu Tyr Lys Cys Lys Val Ser Asn
      195     200     205
30     Lys Ala Leu Pro Ala Pro Ile Glu Lys Thr Ile Ser Lys Ala Lys Gly
      210     215     220
      Gln Pro Arg Glu Pro Gln Val Tyr Thr Leu Pro Pro Ser Arg Glu Glu
      225     230     235     240
35     Met Thr Lys Asn Gln Val Ser Leu Thr Cys Leu Val Lys Gly Phe Tyr
      245     250     255
      Pro Ser Asp Ile Ala Val Glu Trp Glu Ser Asn Gly Gln Pro Glu Asn
      260     265     270
      Asn Tyr Lys Thr Thr Pro Pro Val Leu Asp Ser Asp Gly Ser Phe Phe
      275     280     285
40     Leu Tyr Ser Lys Leu Thr Val Asp Lys Ser Arg Trp Gln Gln Gly Asn
      290     295     300
      Val Phe Ser Cys Ser Val Met His Glu Ala Leu His Asn His Tyr Thr
      305     310     315     320
      Gln Lys Ser Leu Ser Leu Ser Pro Gly Lys
      325     330
45

```

&lt;210&gt; 102

&lt;211&gt; 330

&lt;212&gt; PRT

50 &lt;213&gt; Artificial Sequence

&lt;220&gt;

&lt;223&gt; Human IgG1 heavy chain constant region of IgG1 G1m3 allotype

55 &lt;400&gt; 102

Ala Ser Thr Lys Gly Pro Ser Val Phe Pro Leu Ala Pro Ser Ser Lys  
 1 5 10 15  
 Ser Thr Ser Gly Gly Thr Ala Ala Leu Gly Cys Leu Val Lys Asp Tyr  
 20 25 30  
 Phe Pro Glu Pro Val Thr Val Ser Trp Asn Ser Gly Ala Leu Thr Ser  
 35 40 45  
 Gly Val His Thr Phe Pro Ala Val Leu Gln Ser Ser Gly Leu Tyr Ser  
 50 55 60  
 Leu Ser Ser Val Val Thr Val Pro Ser Ser Ser Leu Gly Thr Gln Thr  
 65 70 75 80  
 Tyr Ile Cys Asn Val Asn His Lys Pro Ser Asn Thr Lys Val Asp Lys  
 85 90 95  
 Arg Val Glu Pro Lys Ser Cys Asp Lys Thr His Thr Cys Pro Pro Cys  
 100 105 110  
 Pro Ala Pro Glu Ala Ala Gly Gly Pro Ser Val Phe Leu Phe Pro Pro  
 115 120 125  
 Lys Pro Lys Asp Thr Leu Met Ile Ser Arg Thr Pro Glu Val Thr Cys  
 130 135 140  
 Val Val Val Asp Val Ser His Glu Asp Pro Glu Val Lys Phe Asn Trp  
 145 150 155 160  
 Tyr Val Asp Gly Val Glu Val His Asn Ala Lys Thr Lys Pro Arg Glu  
 165 170 175  
 Glu Gln Tyr Asn Ser Thr Tyr Arg Val Val Ser Val Leu Thr Val Leu  
 180 185 190  
 His Gln Asp Trp Leu Asn Gly Lys Glu Tyr Lys Cys Lys Val Ser Asn  
 195 200 205  
 Lys Ala Leu Pro Ala Pro Ile Glu Lys Thr Ile Ser Lys Ala Lys Gly  
 210 215 220  
 Gln Pro Arg Glu Pro Gln Val Tyr Thr Leu Pro Pro Ser Arg Glu Glu  
 225 230 235 240  
 Met Thr Lys Asn Gln Val Ser Leu Thr Cys Leu Val Lys Gly Phe Tyr  
 245 250 255  
 Pro Ser Asp Ile Ala Val Glu Trp Glu Ser Asn Gly Gln Pro Glu Asn  
 260 265 270  
 Asn Tyr Lys Thr Thr Pro Pro Val Leu Asp Ser Asp Gly Ser Phe Phe  
 275 280 285  
 Leu Tyr Ser Lys Leu Thr Val Asp Lys Ser Arg Trp Gln Gln Gly Asn  
 290 295 300  
 Val Phe Ser Cys Ser Val Met His Glu Ala Leu His Asn His Tyr Thr  
 305 310 315 320  
 Gln Lys Ser Leu Ser Leu Ser Pro Gly Lys  
 325 330

<210> 103  
 <211> 330  
 <212> PRT  
 <213> Artificial Sequence

<220>  
 <223> Human IgG1 heavy chain constant region of IgG1 G1m3 allotype

<400> 103

	Ala	Ser	Thr	Lys	Gly	Pro	Ser	Val	Phe	Pro	Leu	Ala	Pro	Ser	Ser	Lys
	1				5					10					15	
	Ser	Thr	Ser	Gly	Gly	Thr	Ala	Ala	Leu	Gly	Cys	Leu	Val	Lys	Asp	Tyr
				20					25					30		
5	Phe	Pro	Glu	Pro	Val	Thr	Val	Ser	Trp	Asn	Ser	Gly	Ala	Leu	Thr	Ser
			35					40					45			
	Gly	Val	His	Thr	Phe	Pro	Ala	Val	Leu	Gln	Ser	Ser	Gly	Leu	Tyr	Ser
		50					55					60				
	Leu	Ser	Ser	Val	Val	Thr	Val	Pro	Ser	Ser	Ser	Leu	Gly	Thr	Gln	Thr
10	65					70						75				80
	Tyr	Ile	Cys	Asn	Val	Asn	His	Lys	Pro	Ser	Asn	Thr	Lys	Val	Asp	Lys
				85						90					95	
	Arg	Val	Glu	Pro	Lys	Ser	Cys	Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys
			100						105					110		
15	Pro	Ala	Pro	Glu	Leu	Leu	Gly	Gly	Pro	Ser	Val	Phe	Leu	Phe	Pro	Pro
			115				120						125			
	Lys	Pro	Lys	Asp	Thr	Leu	Met	Ile	Ser	Arg	Thr	Pro	Glu	Val	Thr	Cys
		130					135					140				
	Val	Val	Val	Asp	Val	Ser	His	Glu	Asp	Pro	Glu	Val	Lys	Phe	Asn	Trp
	145					150					155					160
20	Tyr	Val	Asp	Gly	Val	Glu	Val	His	Asn	Ala	Lys	Thr	Lys	Pro	Arg	Glu
				165						170					175	
	Glu	Gln	Tyr	Asn	Ser	Thr	Tyr	Arg	Val	Val	Ser	Val	Leu	Thr	Val	Leu
			180						185					190		
	His	Gln	Asp	Trp	Leu	Asn	Gly	Lys	Glu	Tyr	Lys	Cys	Lys	Val	Ser	Asn
25			195				200						205			
	Lys	Ala	Leu	Pro	Ala	Pro	Ile	Glu	Lys	Thr	Ile	Ser	Lys	Ala	Lys	Gly
		210					215					220				
	Gln	Pro	Arg	Glu	Pro	Gln	Val	Tyr	Thr	Leu	Pro	Pro	Ser	Arg	Glu	Glu
	225					230					235					240
30	Met	Thr	Lys	Asn	Gln	Val	Ser	Leu	Thr	Cys	Leu	Val	Lys	Gly	Phe	Tyr
				245						250					255	
	Pro	Ser	Asp	Ile	Ala	Val	Glu	Trp	Glu	Ser	Asn	Gly	Gln	Pro	Glu	Asn
				260					265						270	
35	Asn	Tyr	Lys	Thr	Thr	Pro	Pro	Val	Leu	Asp	Ser	Asp	Gly	Ser	Phe	Phe
			275					280					285			
	Leu	Tyr	Ser	Lys	Leu	Thr	Val	Asp	Lys	Ser	Arg	Trp	Gln	Gln	Gly	Asn
		290					295					300				
	Val	Phe	Ser	Cys	Ser	Val	Met	His	Glu	Ala	Leu	His	Asn	His	Tyr	Thr
	305					310					315					320
40	Gln	Lys	Ser	Leu	Ser	Leu	Ser	Pro	Gly	Lys						
				325						330						

<210> 104

<211> 107

45 <212> PRT

<213> Artificial Sequence

<220>

<223> Human kappa light chain constant region

50

<400> 104

55

```

Arg Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu
 1          5          10          15
Gln Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe
 20          25          30
5 Tyr Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln
 35          40          45
Ser Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser
 50          55          60
10 Thr Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu
 65          70          75          80
Lys His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser
 85          90          95
Pro Val Thr Lys Ser Phe Asn Arg Gly Glu Cys
100          105

```

15

<210> 105  
 <211> 106  
 <212> PRT  
 <213> Artificial Sequence

20

<220>  
 <223> Human kappa light chain constant region

25

<400> 105

```

Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp Glu Gln
 1          5          10          15
Leu Lys Ser Gly Thr Ala Ser Val Val Cys Leu Leu Asn Asn Phe Tyr
 20          25          30
30 Pro Arg Glu Ala Lys Val Gln Trp Lys Val Asp Asn Ala Leu Gln Ser
 35          40          45
Gly Asn Ser Gln Glu Ser Val Thr Glu Gln Asp Ser Lys Asp Ser Thr
 50          55          60
35 Tyr Ser Leu Ser Ser Thr Leu Thr Leu Ser Lys Ala Asp Tyr Glu Lys
 65          70          75          80
His Lys Val Tyr Ala Cys Glu Val Thr His Gln Gly Leu Ser Ser Pro
 85          90          95
Val Thr Lys Ser Phe Asn Arg Gly Glu Cys
100          105

```

40

<210> 106  
 <211> 990  
 <212> DNA  
 <213> Artificial Sequence

45

<220>  
 <223> Human IgG1 heavy chain constant region of IgG1 G1m3 allotype

50

<400> 106

55

```

5      ggcctccacca aggggtccatc ggtctttcccc ctggcaccctt cctccaagag caacctctggg 60
      ggcacagcgg ccttgggctg cctggtcaag gactacttcc cccaaccggg gacggtgtcg 120
      tggaaactcag gcgccctgac cagcggcgtg cacacettcc cggtctgtcc acagtcctca 180
      ggactctact ccttcagcag cgtggtgacc gtgcctcca gcagcttggg caccacagacc 240
      tacatctgca acgtgaatca caagcccagc aacaccaagg tggacaagag agttgagccc 300
      aaatctttgtg acaaaaactca cacatgcccc ccgtgcccag cacctgaact cctggggggg 360
      ccgtcagttt tctctttccc cccaaaaacc aaggacaccc tcatgatctc cgggaccctt 420
      gaggtcacat gcgtgggtgt ggacgtgagc cacggaagac ctgaggtcaa gttcaactgg 480
      tacgtggacg gcgtggagggt gcataatgcc aagacaaagc cgggggagga gcagtacaac 540
10     agcacgtacc gtgtgggtcag cgtcctcacc gtctctgacc aggactggct gaatggcaag 600
      gagtacaagt gcaagggtct caacaaagcc ctcccagccc ccatogagaa aaccatctcc 660
      aaagccaaag ggcagccccg agaaccacag gtgtacaccc tgcctccacc ccgggaggag 720
      atgaccaaag accaggtcag cctgacctgc ctggtcaaag gcttctatcc cagcgacatc 780
      gccgtggagt gggagagcaa tgggcagccg gagaacaact adaagaccac gctcccctgt 840
15     ctggactcgg accggtcctt cttcctctat agcaagctca ccgtggacaa gagcagggtg 900
      cagcagggga acgtctttct atggtccgtg atgcatgagg ctctgcacaa ccactacacg 960
      cagaagagcc tctccctgtc cccgggtaaa

```

```

20     <210> 107
      <211> 321
      <212> DNA
      <213> Artificial Sequence

      <220>
25     <223> Human kappa light chain constant region

```

```

      <400> 107

30     cgaactgtgg ctgcaccatc tgtcttcacc ttcccggcct ctgatgagca gttgaaatct 60
      ggaactgcct ctgttgtgtg cctgctgaat aacttctatc ccagagaggg caaagtacag 120
      tggaaaggtg ataacgcctt ccaatcgggt aactcccagg agagtgtcac agagcaggac 180
      agcaaggaca gacccctacg cctcagcagc accctgagcg tgagcaaacg agactacgag 240
      aaacacaaaag tctacgcctg cgaagtcccc catcagggcc tgagctcgcc cgtcacaaaag 300
35     agcttcaaca ggggagagtg t

```

```

      <210> 108
      <211> 318
      <212> DNA
      <213> Artificial Sequence

40     <220>
      <223> Human kappa light chain constant region

```

```

45     <400> 108

      actgtggctg caccatctgt ettcattctt ccgccatctg atgagcagtt gaaatctgga 60
      actgcctctg ttgtgtgcct gctgaataac ttctatccc gagaggccaa agtacagtgg 120
      aaggctggatc acgcccctca atcgggtaac tcccaggaga gtgtcacaga gcaggacagc 180
50     aaggacagca cctacagcct cagcagcacc ctgacgctga gcaaaagcaga ctaccagaaa 240
      cacaaagtct acgcctgcga agtcacccat cagggcctga gctcggccct cacaaagagc 300
      ttcacacagg gagagttg

```

```

55     <210> 109
      <211> 147
      <212> PRT
      <213> Homo sapiens

      <400> 109

```

Met Ala Ser His Arg Leu Leu Leu Leu Cys Leu Ala Gly Leu Val Phe  
 1 5 10 15  
 Val Ser Glu Ala Gly Pro Thr Gly Thr Gly Glu Ser Lys Cys Pro Leu  
 20 25 30  
 5 Met Val Lys Val Leu Asp Ala Val Arg Gly Ser Pro Ala Ile Asn Val  
 35 40 45  
 Ala Val His Val Phe Arg Lys Ala Ala Asp Asp Thr Trp Glu Pro Phe  
 50 55 60  
 10 Ala Ser Gly Lys Thr Ser Glu Ser Gly Glu Leu His Gly Leu Thr Thr  
 65 70 75 80  
 Glu Glu Glu Phe Val Glu Gly Ile Tyr Lys Val Glu Ile Asp Thr Lys  
 85 90 95  
 Ser Tyr Trp Lys Ala Leu Gly Ile Ser Pro Phe His Glu His Ala Glu  
 100 105 110  
 15 Val Val Phe Thr Ala Asn Asp Ser Gly Pro Arg Arg Tyr Thr Ile Ala  
 115 120 125  
 Ala Leu Leu Ser Pro Tyr Ser Tyr Ser Thr Thr Ala Val Val Thr Asn  
 130 135 140  
 Pro Lys Glu  
 145

20  
 <210> 110  
 <211> 127  
 <212> PRT  
 <213> Homo sapiens  
 25  
 <400> 110

Gly Pro Thr Gly Thr Gly Glu Ser Lys Cys Pro Leu Met Val Lys Val  
 1 5 10 15  
 Leu Asp Ala Val Arg Gly Ser Pro Ala Ile Asn Val Ala Val His Val  
 20 25 30  
 Phe Arg Lys Ala Ala Asp Asp Thr Trp Glu Pro Phe Ala Ser Gly Lys  
 35 40 45  
 Thr Ser Glu Ser Gly Glu Leu His Gly Leu Thr Thr Glu Glu Glu Phe  
 50 55 60  
 35 Val Glu Gly Ile Tyr Lys Val Glu Ile Asp Thr Lys Ser Tyr Trp Lys  
 65 70 75 80  
 Ala Leu Gly Ile Ser Pro Phe His Glu His Ala Glu Val Val Phe Thr  
 85 90 95  
 40 Ala Asn Asp Ser Gly Pro Arg Arg Tyr Thr Ile Ala Ala Leu Leu Ser  
 100 105 110  
 Pro Tyr Ser Tyr Ser Thr Thr Ala Val Val Thr Asn Pro Lys Glu  
 115 120 125

45  
 <210> 111  
 <211> 127  
 <212> PRT  
 <213> Homo sapiens  
 50  
 <400> 111

55

Gly Pro Thr Gly Thr Gly Glu Ser Lys Cys Pro Leu Met Val Lys Val  
 1 5 10 15  
 Leu Asp Ala Val Arg Gly Ser Pro Ala Ile Asn Val Ala Val His Val  
 20 25 30  
 5 Phe Arg Lys Ala Ala Asp Asp Thr Trp Glu Pro Phe Ala Ser Gly Lys  
 35 40 45  
 Thr Ser Glu Ser Gly Glu Leu His Gly Leu Thr Thr Glu Glu Gln Phe  
 50 55 60  
 10 Val Glu Gly Ile Tyr Lys Val Glu Ile Asp Thr Lys Ser Tyr Trp Lys  
 65 70 75 80  
 Ala Leu Gly Ile Ser Pro Phe His Glu His Ala Glu Val Val Phe Thr

15 Ala Asn Asp Ser Gly Pro Arg Arg Tyr Thr Ile Ala Ala Leu Leu Ser  
 100 105 110  
 Pro Tyr Ser Tyr Ser Thr Thr Ala Val Val Thr Asn Pro Lys Glu  
 115 120 125  
 <210> 112  
 <211> 138  
 <212> PRT  
 <213> Homo sapiens

25 <400> 112  
 Met Ala Ser His Arg Leu Leu Leu Leu Cys Leu Ala Gly Leu Val Phe  
 1 5 10 15  
 Val Ser Glu Ala Gly Pro Thr Gly Thr Gly Glu Ser Lys Cys Pro Leu  
 20 25 30  
 30 Met Val Lys Val Leu Asp Ala Val Arg Gly Ser Pro Ala Ile Asn Val  
 35 40 45  
 Ala Val His Val Phe Arg Lys Ala Ala Asp Asp Thr Trp Glu Pro Phe  
 50 55 60  
 35 Ala Ser Gly Lys Thr Ser Glu Ser Gly Glu Leu His Gly Leu Thr Thr  
 65 70 75 80  
 Glu Glu Glu Phe Val Glu Gly Ile Tyr Lys Val Glu Ile Asp Thr Lys  
 85 90 95  
 Ser Tyr Trp Lys Ala Leu Gly Ile Ser Pro Phe His Glu His Ala Glu  
 100 105 110  
 40 Val Val Phe Thr Ala Asn Asp Ser Gly Pro Arg Arg Tyr Ser Tyr Ser  
 115 120 125  
 Thr Thr Ala Val Val Thr Asn Pro Lys Glu  
 130 135

45 <210> 113  
 <211> 9  
 <212> PRT  
 <213> Artificial Sequence

50 <220>  
 <223> Transthyretin residues 89-97

<400> 113  
 Glu His Ala Glu Val Val Phe Thr Ala  
 1 5

55 <210> 114  
 <211> 15

<212> PRT  
 <213> Artificial Sequence

<220>  
 5 <223> Potential transthyretin immunogen

<400> 114

10 Gly Gly Glu His Ala Glu Val Val Phe Thr Ala Gly Gly Lys Gly  
           1                  5                          10                          15

<210> 115  
 <211> 12  
 <212> PRT  
 15 <213> Artificial Sequence

<220>  
 <223> Potential transthyretin immunogen

20 <400> 115

          Cys Gly Gly Glu His Ala Glu Val Val Phe Thr Ala  
           1                  5                          10

25 <210> 116  
 <211> 12  
 <212> PRT  
 <213> Artificial Sequence

30 <220>  
 <223> Potential transthyretin immunogen

<400> 116

35 Glu His Ala Glu Val Val Phe Thr Ala Cys Gly Gly  
           1                  5                          10

<210> 117  
 <211> 10  
 40 <212> PRT  
 <213> Artificial Sequence

<220>  
 <223> 9D5 CDRH1 - Composite Chothia-Kabat

45 <400> 117

          Gly Phe Thr Phe Ser Ser Tyr Thr Met Ser  
           1                  5                          10

50 <210> 118  
 <211> 10  
 <212> PRT  
 <213> Artificial Sequence

55 <220>  
 <223> 14G8 CDRH1 - Composite Chothia-Kabat

<400> 118

Gly Phe Thr Phe Ser Ser Tyr Thr Met Ser  
1 5 10

5

**Patentkrav**

1. Humaniseret antistof, der binder til transthyretin og omfatter en moden variabel tungkædereion ifølge SEQ ID NO: 11 og en moden variabel letkædereion ifølge SEQ ID NO: 19 eller en moden variabel tungkædereion ifølge SEQ ID NO: 65 og en moden variabel letkædereion ifølge SEQ ID NO: 76.
2. Antistof ifølge krav 1, der udviser en human IgG1-isotype.
3. Antistof ifølge krav 1, der udviser en human IgG2- eller IgG4-isotype.
4. Antistof ifølge krav 1, der er et bindingsfragment, eventuelt hvor bindingsfragmentet er et enkeltkædeantistof-, Fab- eller F(ab')<sub>2</sub>-fragment.
5. Antistof ifølge et hvilket som helst af kravene 1-3, hvor den modne variable letkædereion er kondenseret til en konstant letkædereion, og den modne variable tungkædereion er kondenseret til en konstant tungkædereion.
6. Antistof ifølge krav 5, hvor den konstante tungkædereion er en mutant form af en naturlig human konstant tungkædereion, der har reduceret binding til en Fcγ-receptor i forhold til den naturlige humane konstante tungkædereion.
7. Antistof ifølge krav 5, hvor den konstante tungkædereion er af IgG1-isotypen, eventuelt hvor den modne variable tungkædereion er kondenseret til en konstant tungkædereion med sekvensen ifølge SEQ ID NO: 103 med eller uden det C-terminale lysin, og/eller den modne variable letkædereion er kondenseret til en konstant letkædereion med sekvensen ifølge SEQ ID NO: 104 eller 105.
8. Antistof ifølge krav 1, der omfatter en moden tungkæde med sekvensen ifølge SEQ ID NO: 82 med eller uden det C-terminale lysin og en moden letkæde med sekvensen ifølge SEQ ID NO: 86.
9. Farmaceutisk sammensætning, der omfatter antistoffet ifølge et hvilket som helst foregående krav og en farmaceutisk acceptabel bærer.

**10.** Nukleinsyre eller nukleinsyrer, der koder for tungkæden og letkæden i et antistof som beskrevet i et hvilket som helst af kravene 1-8, eventuelt med en sekvens, der omfatter en hvilken som helst af SEQ ID NOS: 50, 52, 92 og 96.

**11.** Rekombinant ekspressionsvektor, der omfatter nukleinsyren ifølge krav 10 eller rekombinante ekspressionsvektorer, der omfatter nukleinsyrerne ifølge krav 10.

**12.** Værtscelle, der er transformeret med den eller de rekombinante ekspressionsvektor(-er) ifølge krav 11.

**13.** Fremgangsmåde til fremstilling af et antistof ifølge krav 1, hvilken fremgangsmåde omfatter:

- (a) dyrkning af celler, som er transformeret med nukleinsyrer, der koder for de tunge og lette kæder i antistoffet, således at cellerne secernerer antistoffet; og
- (b) oprensning af antistoffet fra celledyrkningsmedier.

**14.** Antistof ifølge et hvilket som helst af kravene 1-8 til anvendelse ved behandling eller frembringelse af profylakse af en transthyretinmedieret amyloidose eller forsinkelse af indtræden af en transthyretinmedieret amyloidose hos et individ.

**15.** In vitro-fremgangsmåde til diagnosticering af en transthyretinmedieret amyloidose hos et individ, hvilken fremgangsmåde omfatter, at en biologisk prøve fra individet bringes i kontakt med en effektiv mængde af antistoffet ifølge et hvilket som helst af kravene 1-8, eventuelt hvor den biologiske prøve og/eller en kontrolprøve er blod, serum, plasma eller fastvæv; eventuelt hvor fastvævet er fra hjertet, det perifere nervesystem, det autonome nervesystem, nyrerne, øjnene eller mave-tarm-kanalen.

# DRAWINGS

	10	20	30	40
m9D5VH	EVKLVESGGGLVQPGGSLKLSCAASGFTFS	YITMSWVRQTPE	42	
Mouse Model 1SEQ_H	EVKLVESGGGLVQPGGSLKLSCAASGFTFS	YITMSWVRQTPE	42	
Mouse Model 1MOK_H	EVKLVESGGGLVQPGGSLKLSCAASGFTFS	YITMSWVRQTPE	42	
Human Acceptor BAC02114	EVQLVESGGGLVQPGGSLRLSCAASGFTFS	SYEMNWVROAPG	42	
Human Acceptor AAX82494	QVQLQESGGGLVQPGGSLKLSCAASGFTFS	SYIGMSWVRQTPD	42	
Hu9D5VHV1	EVQLVESGGGLVQPGGSLRLSCAASGFTFS	SYITMSWVRQAPG	42	
Hu9D5VHV2	EVQLVESGGGLVQPGGSLRLSCAASGFTFS	SYITMSWVRQAPG	42	
Hu9D5VHV2b	EVQLVESGGGLVQPGGSLKLSCAASGFTFS	SYITMSWVRQTPE	42	
Hu9D5VHV3	EVQLVESGGGLVQPGGSLRLSCAASGFTFS	SYITMSWVRQAPG	42	
Hu9D5VHV3b	EVQLVESGGGLVQPGGSLKLSCAASGFTFS	SYITMSWVRQTPG	42	
Hu9D5VHV4	EVQLVESGGGLVQPGGSLRLSCAASGFTFS	SYITMSWVRQAPG	42	
Hu9D5VHV4b	EVQLVESGGGLVQPGGSLKLSCAASGFTFS	SYITMSWVRQAPG	42	
Hu9D5VHV5	EVQLVESGGGLVQPGGSLRLSCAASGFTFS	SYITMSWVRQTPE	42	
	50	60	70	80
m9D5VH	KRLELVAEISNSGDTYYPD	TVKGRFTFSRDNAKNTLYLQMS	84	
Mouse Model 1SEQ_H	KKLEWVAIISKGGSTYYPD	TVKGRFTISRDNKNTLYLQMS	84	
Mouse Model 1MOK_H	KRLEWVAISINNGGRTYYPD	TVKGRFTISRDNKNTLYLQMS	84	
Human Acceptor BAC02114	KGLEWVSIITSSGSTIYYADSVKGRFTISRDN	AKNSLYLQMN	84	
Human Acceptor AAX82494	KRLEWVAITISSGGSTYYPDSVKGRFTISRDN	AKNTLYLQMS	84	
Hu9D5VHV1	KRLELVAEISNSGDTYYPD	TVKGRFTFSRDNAKNSLYLQSN	84	
Hu9D5VHV2	KGLELVAEISNSGDTYYPD	TVKGRFTFSRDNAKNSLYLQSN	84	
Hu9D5VHV2b	KRLELVAEISNSGDTYYPD	TVKGRFTISRDNKNTLYLQMS	84	
Hu9D5VHV3	KGLEWVSEISNSGDTYYPD	TVKGRFTFSRDNAKNSLYLQSN	84	
Hu9D5VHV3b	KRLELVAEISNSGDTYYPD	TVKGRFTISRDNKNTLYLQMS	84	
Hu9D5VHV4	KGLEWVSEISNSGDTYYPD	TVKGRFTISRDNKNSLYLQSN	84	
Hu9D5VHV4b	KRLELVAEISNSGDTYYPD	TVKGRFTISRDNKNTLYLQMS	84	
Hu9D5VHV5	KRLELVAEISNSGDTYYPD	TVKGRFTISRDNKNTLYLQMN	84	

FIG. 1A.1

	90	100	110	120	
m9D5VH	SLKSEDTAMYYCARHYYGGYGGWFFDVWGQTIVTSS				124 (SEQ ID NO.1)
Mouse Model 1SEQ_H	SLKSEDTALYYCAR-GAMFGNDFKYPMDRHWGQTSIVTSS				123 (SEQ ID NO.2)
Mouse Model 1MQK_H	SLKSEDTAMYYCVRHEY-----YAMDYWGQGTIVTSS				118 (SEQ ID NO.62)
Human Acceptor BAC02114	SLRAEDTAVYYCARSGGSRYYIYGMVWGQGTIVTSS				124 (SEQ ID NO.3)
Human Acceptor AAX82494	SLKSEDTAMYYCARLYGYRYF-----DYWGQGTIVTSS				120 (SEQ ID NO.4)
Hu9D5VHV1	SLKAEDTAVYYCARHYYGGYGGWFFDVWGQGTIVTSS				124 (SEQ ID NO.5)
Hu9D5VHV2	LLRAEDTAVYYCARHYYGGYGGWFFDVWGQGTIVTSS				124 (SEQ ID NO.6)
Hu9D5VHV2b	SLKSEDTAMYYCARHYYGGYGGWFFDVWGQGLIVTSS				124 (SEQ ID NO.7)
Hu9D5VHV3	LLRAEDTAVYYCARHYYGGYGGWFFDVWGQGLIVTSS				124 (SEQ ID NO.8)
Hu9D5VHV3b	SLKSEDTAMYYCARHYYGGYGGWFFDVWGQGLIVTSS				124 (SEQ ID NO.9)
Hu9D5VHV4	LLRAEDTAVYYCARHYYGGYGGWFFDVWGQGTIVTSS				124 (SEQ ID NO.10)
Hu9D5VHV4b	SLKSEDTAMYYCARHYYGGYGGWFFDVWGQGLIVTSS				124 (SEQ ID NO.11)
Hu9D5VHV5	LLRAEDTAVYYCARHYYGGYGGWFFDVWGQGTIVTSS				124 (SEQ ID NO.12)

FIG. 1A.2

m9D5VL	10	20	30	40
Mouse Model 1MJU L	DIVMTQAAPSPVPPGESVSI			
Human Acceptor ABC66952	DIVMTQAAPSPVPPGESVSI			
Hu9D5VLv1	DIVMTQSPLSLPTPGEPA			
Hu9D5VLv2	DIVMTQSPLSLPTPGEPA			
Hu9D5VLv3	DIVMTQSPLSLPTPGEPA			
Hu9D5VLv4	DIVMTQSPLSLPTPGEPA			
Hu9D5VLv5	DIVMTQSPLSLPTPGEPA			
m9D5VL	50	60	70	80
Mouse Model 1MJU L	FLQKPGQSPQLLIYRVSN			
Human Acceptor ABC66952	FLQKPGQSPQLLIYRVSN			
Hu9D5VLv1	FLQKPGQSPQLLIYRVSN			
Hu9D5VLv2	FLQKPGQSPQLLIYRVSN			
Hu9D5VLv3	FLQKPGQSPQLLIYRVSN			
Hu9D5VLv4	FLQKPGQSPQLLIYRVSN			
Hu9D5VLv5	FLQKPGQSPQLLIYRVSN			
m9D5VL	90	100	110	
Mouse Model 1MJU L	SRVEAEDVGVVYVCMQH			
Human Acceptor ABC66952	SRVEAEDVGVVYVCMQH			
Hu9D5VLv1	SRVEAEDVGVVYVCMQH			
Hu9D5VLv2	SRVEAEDVGVVYVCMQH			
Hu9D5VLv3	SRVEAEDVGVVYVCMQH			
Hu9D5VLv4	SRVEAEDVGVVYVCMQH			
Hu9D5VLv5	SRVEAEDVGVVYVCMQH			

FIG. 1B

(SEQ ID NO.16)  
 (SEQ ID NO.17)  
 (SEQ ID NO.18)  
 (SEQ ID NO.19)  
 (SEQ ID NO.20)  
 (SEQ ID NO.21)  
 (SEQ ID NO.22)  
 (SEQ ID NO.23)






FIG. 2A

	10	20	30	40
m14G8VL	DIVMTQAAPSPVPTPGESVVISCRSNKSKLLHSHNGNTYLYW 40			
Mouse Model 1MJU_L	DIVMTQAAPSPVPTPGESVVISCRSNKSKLLHSHNGNTYLYW 40			
Human Acceptor ABA71374	DIVMTQTPPLSLPVTGEPASISCRSSQSLRHSHYGYIYDW 40			
Human Acceptor ABC66952	DIVMTQSPPLSLPVTGEPASISCRSNOSLHSHNGYINLWDW 40			
Hu14G8VLv1	DIVMTQAPSLPVTGEPESVVISCRSNKSKLLHSHNGNTYLYW 40			
Hu14G8VLv2	DIVMTQSPPLSLPVTGEPASISCRSNKSKLLHSHNGNTYLYW 40			
Hu14G8VLv3	DIVMTQSPPLSLPVTGEPASISCRSSKSKLLHSHNGNTYLYW 40			
	50	60	70	80
m14G8VL	FLQRPQSPQLLIYRVSNNLASGVPPDRFSGSGGTAFTLRI 80			
Mouse Model 1MJU_L	FLQRPQSPQLLIYRMSNNLASGVPPDRFSGSGGTAFTLRI 80			
Human Acceptor ABA71374	YLQKPGQSPQVLIYLGSNRASGVPPDRFSGSGGTDFTLKI 80			
Human Acceptor ABC66952	YLQKPGQSPQLLIYSGSNRASGVPPDRFSGSGGTDFTLKI 80			
Hu14G8VLv1	FLQKPGQSPQLLIYRVSNNLASGVPPDRFSGSGGTAFTLKI 80			
Hu14G8VLv2	FLQKPGQSPQLLIYRVSNNLASGVPPDRFSGSGGTDFTLKI 80			
Hu14G8VLv3	FLQKPGQSPQLLIYRVSNNLASGVPSRFSGSGGTDFTLKI 80			
	90	100	110	
m14G8VL	SRVEAEDGVVYVCMQHLEYPLTFGAGTKLELKR 113 (SEQ ID NO.70)			
Mouse Model 1MJU_L	SRVEAEDGVVYVCMQHLEYPLTFGAGTKLELKR 113 (SEQ ID NO.71)			
Human Acceptor ABA71374	SRVEAEDGVVYVCMALQVCSFGQGTKEIKR 113 (SEQ ID NO.72)			
Human Acceptor ABC66952	SRVEAEDGVVYVCMALQSPYTFGQGTKEIKR 113 (SEQ ID NO.73)			
Hu14G8VLv1	SRVEAEDGVVYVCMQHLEYPLTFGQGTKEIKR 113 (SEQ ID NO.74)			
Hu14G8VLv2	SRVEAEDGVVYVCMQHLEYPLTFGQGTKEIKR 113 (SEQ ID NO.75)			
Hu14G8VLv3	SRVEAEDGVVYVCMQHLEYPLTFGQGTKEIKR 113 (SEQ ID NO.76)			

FIG. 2B

**FIG. 3A****FIG. 3B**




FIG. 4A




FIG. 4B




FIG. 4C




FIG. 5A




FIG. 5B




FIG. 6




FIG. 7A




FIG. 7B