

(19) **DANMARK**

(10) **DK/EP 2967836 T3**



(12) **Oversættelse af  
europæisk patentskrift**

Patent- og  
Varemærkestyrelsen

- 
- (51) Int.Cl.: **A 61 F 2/10 (2006.01)** **A 01 N 1/02 (2006.01)** **A 61 F 2/00 (2006.01)**
- (45) Oversættelsen bekendtgjort den: **2019-01-28**
- (80) Dato for Den Europæiske Patentmyndigheds bekendtgørelse om meddelelse af patentet: **2018-10-10**
- (86) Europæisk ansøgning nr.: **14775892.4**
- (86) Europæisk indleveringsdag: **2014-03-13**
- (87) Den europæiske ansøgnings publiceringsdag: **2016-01-20**
- (86) International ansøgning nr.: **US2014025599**
- (87) Internationalt publikationsnr.: **WO2014160000**
- (30) Prioritet: **2013-03-13 US 201361779661 P**
- (84) Designerede stater: **AL AT BE BG CH CY CZ DE DK EE ES FI FR GB GR HR HU IE IS IT LI LT LU LV MC MK MT NL NO PL PT RO RS SE SI SK SM TR**
- (73) Patenthaver: **Stratatech Corporation, 505 South Rosa Road Suite 169, Madison, Wisconsin 53719, USA**
- (72) Opfinder: **ALLEN-HOFFMANN, B. Lynn, 3905 Council Crest, Madison, Wisconsin 53711, USA**  
**PIRNSTILL, John C., 2765 Jacquelyn Drive, Fitchburg, Wisconsin 53711, USA**  
**GRATZ, Kenneth R., 733 Miami Pass, Madison, Wisconsin 53711, USA**  
**COMER, Allen R., 5809 Chester Circle, Madison, Wisconsin 53719, USA**
- (74) Fuldmægtig i Danmark: **Zacco Denmark A/S, Arne Jacobsens Allé 15, 2300 København S, Danmark**
- (54) Benævnelse: **Kryokonservering af levedygtige menneskehudssubstitutter**
- (56) Fremdragne publikationer:  
**WO-A1-95/07611**  
**WO-A1-96/24018**  
**WO-A1-96/33750**  
**WO-A1-97/28402**  
**WO-A1-2008/088136**  
**WO-A1-2009/065005**  
**WO-A2-96/14738**  
**WO-A2-02/058588**  
**US-A- 5 891 617**  
**US-A- 6 110 208**  
**US-A1- 2004 053 409**  
**US-A1- 2010 119 615**  
**US-A1- 2013 108 670**  
**US-A1- 2013 108 683**  
**HANS-JÜRGEN STARK ET AL: "Organotypic cocultures as skin equivalents: A complex and sophisticated in vitro system", BIOLOGICAL PROCEDURES ONLINE, vol. 19, no. 4, 12 April 2004 (2004-04-12) , pages 55-60,**

Fortsættes ...

XP055268292, DOI: 10.1251/bpo72

T Rupf ET AL: "Cryopreservation of Organotypical Cultures Based on 3D Scaffolds", *Cryoletters*, 1 March 2010 (2010-03-01), pages 157-168, XP055285494, England Retrieved from the Internet:

URL:[http://www.ingentaconnect.com/content/cryo/cryo/2010/00000031/00000002/art00008?](http://www.ingentaconnect.com/content/cryo/cryo/2010/00000031/00000002/art00008?token=004b1623dac778437a63736a6f5e47745d36664670502e796f644a467b4d616d3f4e4b34a25)

token=004b1623dac778437a63736a6f5e47745d36664670502e796f644a467b4d616d3f4e4b34a25

C C BONDOC ET AL: "Clinical Experience with Viable Frozen Human Skin and a Frozen Skin Bank", *ANN. SURG.*, vol. 174, no. 3, 1 September 1971 (1971-09-01), pages 371-381, XP055285286,

J Demetrulias ET AL: "Skin2(R) - an in vitro human skin model: the correlation between in vivo and in vitro testing of surfactants", *Exp Dermatol*, 1 January 1998 (1998-01-01), pages 18-26, XP055286322, Retrieved from the Internet: URL:<http://onlinelibrary.wiley.com/store/10.1111/j.1600-0625.1998.tb00298.x/asset/j.1600-0625.1998.tb00298.x.pdf?v=1&t=iqas0eww&s=5689b9b94b21794735fa1894ac95049a933f5727> [retrieved on 2016-07-06]

BONDOC, CC ET AL.: 'Clinical Experience With Viable Frozen Human Skin And A Frozen Skin Bank.' *ANN. SURG.* vol. 174, no. 3, 1971, pages 371 - 381, XP055285286

# DESCRIPTION

## FIELD OF THE INVENTION

**[0001]** The present invention relates generally to systems and methods for cryopreservation of viable human skin substitutes.

## BACKGROUND

**[0002]** The emerging field of tissue engineering (TE) is poised to make enormous progress in the treatment of organ disease and dysfunction in the coming decade. In 2001, there were 23 cell-based therapeutics approved for market in the United States (U.S.) and Europe, of which nine were skin substitutes or grafts, and 100 more products were in development. (De Bree, Genomics-based Drug Data Report and Regenerative Therapy (1) 2:77-96 (2001)). A decade later, engineered tissues have emerged as a discrete industry sector within the wound care industry and represent one of regenerative medicine's most promising cell based therapeutic platforms. The global wound care market was estimated to be \$16.8 billion in 2012 and has been growing at a rate of approximately 6% annually (Kalorama Information, April 2012). The bulk of this market is comprised of traditional sectors that are mature and highly competitive, and include products targeting basic and advanced wound care, wound closure, and anti-infectives. This has led competitors to increasingly focus their attention on developing highly differentiated products in the more innovative active wound care sector, a sector that represents approximately 15% of the overall market. While sales of negative pressure wound therapy still predominate, 2010 sales of engineered tissues and other products within the U.S. biologics market grew to \$448 million, and are projected to increase to \$1.058 billion by 2015, a compound annual growth rate of 18.8% (BioMedGPS-SmarTRAK market analysis, 2012).

**[0003]** Although a multitude of revolutionary and economically important applications for engineered tissues and organs exist in the human health arena, the full economic potential of the industry is far from realized. At present, only two tissue engineering companies worldwide have been able to commercialize cell based, skin substitute products focused on cutaneous wound healing and achieve annual sales in excess of \$100 million.

**[0004]** A major impediment to the acceptance of engineered tissues by medical practitioners, healthcare providers, and second party payers is the lack of a means to effectively and efficiently preserve and store engineered tissues. The nature of living cells and tissue products makes development of long-term storage challenging. Current engineered tissues must often be stored and shipped under carefully controlled conditions to maintain viability and function. Typically, engineered tissue products take weeks or months to produce but must be used within hours or days after manufacture. As a result, TE companies must continually operate with their production facilities at top capacity and absorb the costs of unsold product which

must be discarded. These inventory losses, on top of already costly manufacturing process, have forced prices to impractical levels. As one specific example, APLIGRAF requires about four weeks to manufacture, is usable for only 15 days and must be maintained between 20 and 23°C until used. As another example, EPICEL is transported by a nurse from Genzyme Biosurgery's production facility in Cambridge, MA to the point of use in a portable incubator and is used immediately upon arrival. Such constraints represent significant challenges to developing convenient and cost-effective products.

**[0005]** Cryopreservation has been explored as a solution to the storage problem, but it is known to induce tissue damage through ice formation, chilling injury, and osmotic imbalance. Besides APLIGRAF, the only other approved living skin equivalent, ORCEL, has been evaluated as a frozen product but had the drawback that it must be maintained at temperatures below -100°C prior to use. This requires specialized product delivery and storage conditions, including the use of dangerous goods during transport, and use of liquid nitrogen for storage, which is expensive, dangerous, and not readily available in rural clinics and field hospitals.

**[0006]** Accordingly, what is needed in the art are improved methods of cryopreserving viable engineered tissues and cells for storage under conditions that are routinely available at the point of use.

#### **SUMMARY OF THE INVENTION**

**[0007]** The invention is set out in the appended claims. Namely, the invention is directed to a method of cryopreserving an organotypically cultured skin equivalent to maintain viable tissue comprising: a) treating an organotypically cultured skin equivalent in a cryoprotectant solution in a single step, wherein said organotypically cultured skin equivalent comprises stratified squamous epithelia on a dermal layer comprising fibroblast, and wherein said cryoprotectant is provided in a solution comprising 21% to 70% of said solution by volume and said cryoprotectant is glycerol; b) separating the treated organotypically cultured skin equivalent from excess cryoprotectant solution, and packaging said treated organotypically cultured skin equivalent in the absence of additional cryoprotectant to provide a packaged skin equivalent, wherein said packaging preferably further comprises enclosing said cryopreserved skin equivalent in a sterile bag and enclosing said sterile bag in a second bag; and c) freezing said packaged organotypically cultured skin equivalent to provide a cryopreserved skin equivalent.

**[0008]** The embodiments of the description which do not fall within the scope of said claims are provided for illustrative purposes only and do not form part of the present invention. The present invention relates generally to systems and methods for cryopreservation of viable human skin substitutes.

**[0009]** The present description provides methods of cryopreserving an organotypically cultured skin equivalent to maintain viable tissue comprising: treating an organotypically cultured skin

equivalent in a cryoprotectant solution in a single step; packaging the organotypically cultured skin equivalent to provide a packaged skin equivalent; and freezing the organotypically cultured skin equivalent to provide a packaged cryopreserved skin equivalent. The cryoprotectant may be provided in a solution comprising about 20% or 21% to about 70% of the solution by volume, and more preferably about 20% or 21% to about 45% of the solution by volume or 37.5% to 62.5% of the solution by volume, or most preferably from about 25% to 40% of the solution by volume or 42.5% to 57.5% of the solution by volume, depending on the temperature. The treatment with cryoprotectant may be conducted at from about 2 °C to 8 °C, while in other examples, the treatment step is conducted at room temperature, for example from about 15 °C to 30 °C. According to the invention, the cryoprotectant is glycerol. The freezing may further comprise freezing the organotypically cultured skin equivalent in the absence of substantial excess cryoprotectant. The freezing may further comprise freezing at about -80 °C. The freezing may further comprise direct exposure to temperatures ranging from about -50 °C to -100 °C. The packaging may further comprise enclosing the cryopreserved skin equivalent in a sterile bag and enclosing the sterile bag in a second bag. The organotypically cultured skin equivalents may comprise NIKS cells. The NIKS cells may comprise an exogenous nucleic acid sequence encoding an exogenous polypeptide. The skin equivalent may retain viability after thawing. The skin equivalent may have an  $A_{550}$  of at least 50% of a reference skin equivalent as determined by an MTT assay.

**[0010]** The methods may further comprise thawing said cryopreserved skin equivalent and applying said thawed skin equivalent to a patient in need thereof, wherein said thawed skin equivalent is not rinsed prior to said application to said patient. The present description provides methods of thawing a cryopreserved skin equivalent prior to application to a subject, comprising: warming the cryopreserved skin equivalent; and contacting the cryopreserved skin equivalent with a diffusion mediator comprising a tissue compatible solution to allow removal of the cryoprotectant solution by diffusion. The warming may comprise exposure to room temperature at the site of use. The diffusion mediator may be selected from the group consisting of an absorbent medium, a membrane, and a dialysis bag. The absorbent medium may be selected from the group consisting of Telfa pads, foam pads, gauze pads, and cellulosic pads containing the tissue compatible medium. The tissue compatible solution may be a buffered solution.

**[0011]** The present description provides methods of treating a subject comprising providing a packaged cryopreserved skin equivalent produced as described above; aseptically transferring the cryopreserved skin equivalent from the package; warming the cryopreserved skin equivalent; contacting the cryopreserved skin equivalent with an absorbent medium comprising a tissue compatible solution to allow removal of the cryoprotectant solution by diffusion; and applying the cryopreserved skin equivalent to the subject.

**[0012]** The present description provides a cryopreserved skin equivalent equilibrated with a cryoprotectant, the skin equivalent being substantially free of excess cryoprotectant on the exterior surface of the skin equivalent. The present description provides a system comprising the foregoing skin equivalent disposed on an absorbent medium.

**[0013]** The present description provides methods comprising providing the packaged cryopreserved skin equivalent as described above; and applying the skin equivalent to a wound under conditions such that the skin equivalent contacts the wound.

**[0014]** The present description provides a kit comprising a cryopreserved skin substitute, an absorbent medium, and a tissue compatible solution. The cryopreserved skin substitute may be packaged in a sealable enclosure. The cryopreserved skin substitute may be provided in a culture vessel packaged in the bag.

**[0015]** The present description provides a method comprising: providing a culture dish comprising a cell culture substrate movable between defined upper and lower positions in the culture dish, forming a skin equivalent on the cell culture substrate, wherein the cell culture substrate is at the upper position, lowering the cell culture substrate to the lower position for further processing. The further processing may comprise treating the skin equivalent with a cryoprotectant solution. The further processing may comprise freezing the skin equivalent in the culture dish.

**[0016]** The present description provides a method of producing a cryopreserved skin equivalent comprising: providing a culture dish comprising an insert movable between upper and lower positions in the culture dish, the insert having a bottom planar surface formed from a porous membrane, forming a dermal equivalent comprising fibroblast cells on the porous membrane in the insert, wherein the insert is placed the upper position the culture dish, culturing the fibroblast cells to form a dermal equivalent, applying keratinocyte cells to the dermal equivalent, culturing the keratinocytes in a culture medium under conditions such that the keratinocytes form a skin equivalent comprising stratified epithelium, removing the culture medium, lowering the insert to the lower position, treating the skin equivalent with a cryoprotectant solution, and freezing the skin equivalent in the culture dish.

**[0017]** The present disclosure relates to methods of treating a patient in need thereof with a cryopreserved skin equivalent made by the foregoing methods comprising thawing said cryopreserved skin equivalent and applying said thawed skin equivalent to said patient in need thereof, wherein said thawed skin equivalent is not rinsed prior to said application to said patient. For the sake of clarity, this aspect of the present disclosure is not part of the invention as defined in the appended claims.

## **DESCRIPTION OF FIGURES**

### **[0018]**

FIG. 1 is a graph of viability testing of thawed tissues using the MTT assay. The tissues were cryopreserved using a controlled-rate freezer (CRF) or by passive freezing in an ultra-cold freezer (set to -80 °C) and stored for up to 3 months, then thawed and analyzed at 1 hour and

1 day. Viability was measured by MTT analysis of three 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue at each time point (mean +/- 1 st.dev.).

FIG. 2 is a temperature profile from a simulated passive tissue freeze. Temperature probes were affixed to the bottom of three 100 mm x 20 mm dishes with a Transwell insert on top. The lid was placed on the dish and then they were packaged within an inner Whirl-pak pouch and an outer Mylar pouch and placed in an ultra-cold freezer set to -80 °C. The packaged tissues were placed in a metal rack at a bottom, middle, and top position.

FIG. 3 is a graph of viability testing of thawed tissues using the MTT assay. The tissues were cryopreserved using a CRF and then stored in ultra-cold freezers set to -80 °C or -50 °C for 1 to 6 months, then thawed and analyzed at 1 hour and 1 day. Viability was measured by MTT analysis of three 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue at each time point (mean +/- 1 st.dev.).

FIG. 4 is a graph of viability testing of thawed tissues using the MTT assay. Tissues from two independent lots were cryopreserved and then stored in an ultra-cold freezer set to -80 °C for greater than 1 month. Two tissues from each lot were placed into a container of dry ice for 48 hours, with two control tissues from each lot remaining in the ultra-cold freezer. At the end of the dry ice exposure, all tissues were thawed and analyzed at 15 minutes and 1 day. Viability was measured by MTT analysis of four 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue at each time point (mean +/- 1 st.dev.).

FIG. 5 is a graph of viability testing of thawed tissues using the MTT assay. The tissues were cryopreserved using a CRF following a stepwise temperature reduction (room temperature to 2-8 °C to -20 °C) during which the tissues were exposed to a graded series of glycerol concentrations soaked into sterile cellulose pads. All tissues were initially exposed to 16.3% glycerol at room temperature. Next, all tissues were switched to 2-8 °C, with two tissues transferred to 32.5% glycerol. Lastly, all tissues were switched to -20 °C, with one tissue transferred to 65% glycerol. Tissues were stored in vapor-phase LN<sub>2</sub> for 6 days. Viability was measured by MTT analysis of three 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue at 1 hour and 1 day post-thaw (mean +/- 1 st.dev.).

FIG. 6 is a graph of viability testing of thawed tissues using the MTT assay. The tissues were cryopreserved in 65% glycerol by passive freezing in an ultra-cold freezer set to -80 °C and then stored at -80 °C for 6 weeks. Tissues were thawed and placed in a hold chamber, consisting of two cellulose filter pads on a raised stainless steel lifter, containing 100 ml of the specified solution and held for 1 hour at the specified temperature. Tissues were analyzed after 1 hour and 1 day. Viability was measured by MTT analysis of three 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue at each time point (mean +/- 1 st.dev.).

FIG. 7 is a graph of viability testing of thawed tissues using the MTT assay. The tissues were cryopreserved in 65% glycerol and then stored in vapor-phase LN<sub>2</sub> for 2 weeks. One tissue was thawed directly into a growth chamber without a cellulose pad, while the other was placed in a hold chamber, consisting of 2 cellulose filter pads on a raised stainless steel lifter,

containing 100 ml of the growth media. The thawed tissues were held for 1 hour at 37 °C. Tissues were analyzed after 1 hour and 1 day. Viability was measured by MTT analysis of three 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue at each time point (mean +/- 1 st.dev.).

FIG. 8 is a graph of viability testing of thawed tissues using the MTT assay. Tissues from one lot were cryopreserved and then stored in an ultra-cold freezer set to -80 °C for 1 week. The tissues were thawed for 10 minutes, placed in hold chambers containing either Telfa pads or Whatman pads soaked with 40 ml of a buffered nutrient solution warmed to either 37 °C or room temperature (n=2 tissues per thaw condition). The tissues were placed in the hold chambers for 15-20 minutes, and then re-cultured in a growth chamber for 1 day. The tissues were analyzed at 1 day post-thaw. Viability was measured by MTT analysis of four 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue (mean +/- 1 st.dev.). The viability specification is indicated by the red dashed line.

FIG. 9 is a graph of viability testing of thawed tissues using the MTT assay. Tissues from three independent lots were treated with 50% glycerol at 2 to 8 °C and cryopreserved as described in Example 2. Cryopreserved tissues were stored in an ultra-cold freezer at approximately -80 °C for up to 12 months. After 0, 2, 3, 5, and 12 months of storage, two tissues from each lot were thawed at room temperature and held on media saturated Telfa pads for 15 minutes. Tissues were then returned to culture overnight in a culture dish with a raised stainless steel lifter, containing 90 ml of growth media. Viability was measured by MTT analysis of four 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue at each time point (mean +/- 1 st.dev.).

FIG. 10 is a graph of viability testing of thawed tissues using the MTT assay. Tissues from one lot were treated with the specified glycerol concentration (32.5% or 50%) at the listed conditions prior to the freeze and then stored in an ultra-cold freezer set to -80 °C for 9 days. The tissues were thawed for 10 minutes, placed on a Telfa stack soaked with 40 ml of a buffered nutrient solution for 15-20 minutes, and then re-cultured in a growth chamber for 1 day. The tissues were analyzed at 1 day post-thaw. Viability was measured by MTT analysis of four 0.5 cm<sup>2</sup> punches from each 44 cm<sup>2</sup> tissue (mean +/- 1 st.dev.). The viability specification is indicated by the red dashed line.

## DEFINITIONS

**[0019]** As used herein, the terms "skin equivalent", "human skin equivalent", "human skin substitute", and "organotypic cultures" are used interchangeably to refer to an in vitro derived culture of keratinocytes that has stratified into squamous epithelia. Typically, the skin equivalents are produced by organotypic culture and include a dermal layer in addition to a keratinocyte layer.

**[0020]** As used herein, the term "sterile" refers to a skin equivalent that is essentially or



completely free of detectable microbial or fungal contamination.

**[0021]** As used herein, the term "NIKS cells" refers to cells having the characteristics of the cells deposited as cell line ATCC CRL-1219. NIKS stands for near-diploid immortalized keratinocytes.

**[0022]** As used herein, the term "viable" when used in reference to a skin equivalent refers to the viability of cells in the skin equivalent following cryopreservation. In preferred embodiments, a "viable" skin has an A<sub>550</sub> of at least 50%, 60%, 70%, 80% or 90% of a control non-cryopreserved tissue as measured by an MTT assay or at least 50%, 60%, 70%, 80% or 90% of the readout value of a similar viability assay.

### **DETAILED DESCRIPTION**

**[0023]** The invention is set out in the appended claims. Namely, the invention is directed to a method of cryopreserving an organotypically cultured skin equivalent to maintain viable tissue comprising: a) treating an organotypically cultured skin equivalent in a cryoprotectant solution in a single step, wherein said organotypically cultured skin equivalent comprises stratified squamous epithelia on a dermal layer comprising fibroblast, and wherein said cryoprotectant is provided in a solution comprising 21% to 70% of said solution by volume and said cryoprotectant is glycerol; b) separating the treated organotypically cultured skin equivalent from excess cryoprotectant solution, and packaging said treated organotypically cultured skin equivalent in the absence of additional cryoprotectant to provide a packaged skin equivalent, wherein said packaging preferably further comprises enclosing said cryopreserved skin equivalent in a sterile bag and enclosing said sterile bag in a second bag; and c) freezing said packaged organotypically cultured skin equivalent to provide a cryopreserved skin equivalent.

**[0024]** The embodiments of the description which do not fall within the scope of said claims are provided for illustrative purposes only and do not form part of the present invention. The present invention relates generally to systems and methods for cryopreservation of human skin substitutes. In particular, the present invention relates to methods for cryopreserving viable human skin equivalents so that they can be stored for prolonged periods at the site of use, such as a hospital, operating room or burn unit. The methods disclosed herein allow for novel increases in efficiency of the preservation process and utilization of preserved skin equivalents, including single-step equilibration in a cryoprotectant, packaging of the skin equivalent in a sterile package prior to cryopreservation, and the ability to use the cryopreserved tissues for direct application to a patient (e.g., in a grafting procedure) without rinsing. A ready-to-use cryopreserved skin equivalents for use in treatment of a patient is provided, preferably for use in grafting procedures. The cryopreserved skin equivalent is designed for long term storage at the site of use. The cryopreserved equivalents may be engineered to deliver the broad spectrum human host defense peptides such as  $\beta$ -defensin-3 (hBD-3) or cathelicidin (hCAP18/LL-37), or pro-angiogenic factors, to the wound bed.

**[0025]** Previously, cadaver skin has been harvested and cryopreserved by treatment with from 10% to 20% glycerol as a cryopreservative. See e.g., Kagan et al., Clin Lab Med 25 (2005) 587-605. Surprisingly, it has been found that increased glycerol concentrations are needed to cryopreserve human skin equivalents.

**[0026]** Accordingly, in some examples, the present description provides a cryopreserved skin equivalent. In some examples, the skin equivalent has been engineered to express and provide exogenous antimicrobial polypeptides or pro-angiogenic factors. There is no limitation to the use of any particular antimicrobial polypeptide. The antimicrobial polypeptide may be human  $\beta$ -defensin-1, human  $\beta$ -defensin-2, human  $\beta$ -defensin-3, or cathelicidin (hCAP-18/LL37) or variant. Nucleic acid constructs or vectors encoding the antimicrobial polypeptide or pro-angiogenic factor may be introduced into the keratinocytes (e.g., NIKS cells) and the transfected keratinocytes are used to make the skin equivalent by organotypic culture techniques. Examples for the production of skin equivalents expressing exogenous polypeptides, as well as additional wild-type and variant antimicrobial polypeptides can be found in U.S. Pat. Nos. 7,674,291; 7,807,148; 7,915,042; 7,988,959; and 8,092,531.

**[0027]** The cryopreserved skin equivalents may be applied to wounds after thawing and left in place. The cryopreserved skin equivalents may be applied temporarily to wounds. The cryopreserved skin equivalents may be removed and replaced with additional cryopreserved human skin equivalents.

#### **A) Skin equivalents produced by organotypic culture**

**[0028]** There is not limitation to the use of any particular source of cells that are capable of differentiating into squamous epithelia. Indeed, the use of a variety of cell lines and sources that can differentiate into squamous epithelia, including both primary and immortalized keratinocytes is contemplated. Sources of cells include keratinocytes and dermal fibroblasts biopsied from humans and cadaveric donors (Auger et al., In Vitro Cell. Dev. Biol. - Animal 36:96-103; U.S. Pat. Nos. 5,968,546 and 5,693,332), neonatal foreskins (Asbill et al., Pharm. Research 17(9): 1092-97 (2000); Meana et al., Burns 24:621-30 (1998); U.S. Pat. Nos. 4,485,096; 6,039,760; and 5,536,656), and immortalized keratinocytes cell lines such as NM1 cells (Baden, In Vitro Cell. Dev. Biol. 23(3):205-213 (1987)), HaCaT cells (Boucamp et al., J. cell. Boil. 106:761-771 (1988)); and NIKS cells (Cell line BC-1-Ep/SL; U.S. Pat. No. 5,989,837; ATCC CRL-12191). Each of the mentioned cell lines can be cultured or genetically modified in order to produce a cell line capable of expressing or co-expressing the desired protein(s). NIKS cells may be utilized. The discovery of the novel NIKS human keratinocyte cell line provides an opportunity to genetically engineer human keratinocytes with non-viral vectors. A unique advantage of the NIKS cells is that they are a consistent source of genetically-uniform, pathogen-free human keratinocytes. For this reason, they are useful for the application of genetic engineering and genomic gene expression approaches to provide human skin equivalents with enhanced properties over currently available skin equivalents. NIKS cells,

identified and characterized at the University of Wisconsin, are nontumorigenic, karyotypically stable, and exhibit normal growth and differentiation both in monolayer and organotypic culture. NIKS cells form fully stratified skin equivalents in culture. These cultures are indistinguishable by all criteria tested thus far from organotypic cultures formed from primary human keratinocytes. Unlike primary cells however, NIKS cells exhibit an extended lifespan in monolayer culture. This provides an opportunity to genetically manipulate the cells and isolate new clones of cells with new useful properties (Allen-Hoffmann et al., *J. Invest. Dermatol.*, 114(3): 444-455 (2000)).

**[0029]** The NIKS cells arose from the BC-1- Ep strain of human neonatal foreskin keratinocytes isolated from an apparently normal male infant. In early passages, the BC-1- Ep cells exhibited no morphological or growth characteristics that were atypical for cultured normal human keratinocytes. Cultivated BC-1- Ep cells exhibited stratification as well as features of programmed cell death. To determine replicative lifespan, the BC-1- Ep cells were serially cultivated to senescence in standard keratinocyte growth medium at a density of  $3 \times 10^5$  cells per 100-mm dish and passaged at weekly intervals (approximately a 1:25 split). By passage 15, most keratinocytes in the population appeared senescent as judged by the presence of numerous abortive colonies which exhibited large, flat cells. However, at passage 16, keratinocytes exhibiting a small cell size were evident. By passage 17, only the small-sized keratinocytes were present in the culture and no large, senescent keratinocytes were evident. The resulting population of small keratinocytes that survived this putative crisis period appeared morphologically uniform and produced colonies of keratinocytes exhibiting typical keratinocyte characteristics including cell-cell adhesion and apparent squame production. The keratinocytes that survived senescence were serially cultivated at a density of  $3 \times 10^5$  cells per 100-mm dish. Typically the cultures reached a cell density of approximately  $8 \times 10^6$  cells within 7 days. This stable rate of cell growth was maintained through at least 59 passages, demonstrating that the cells had achieved immortality. The keratinocytes that emerged from the original senescencing population are now termed NIKS. The NIKS cell line has been screened for the presence of proviral DNA sequences for HIV-1, HIV-2, EBV, CMV, HTLV-1, HTLV-2, HBV, HCV, B-19 parvovirus, HPV-16, SV40, HHV-6, HHV-7, HPV-18 and HPV-31 using either PCR or Southern analysis. None of these viruses were detected.

**[0030]** Chromosomal analysis was performed on the parental BC-1- Ep cells at passage 3 and NIKS cells at passages 31 and 54. The parental BC-1- Ep cells have a normal chromosomal complement of 46, XY. At passage 31, all NIKS cells contained 47 chromosomes with an extra isochromosome of the long arm of chromosome 8. No other gross chromosomal abnormalities or marker chromosomes were detected. The karyotype of the NIKS cells has been shown to be stable to at least passage 54.

**[0031]** The DNA fingerprints for the NIKS cell line and the BC-1- Ep keratinocytes are identical at all twelve loci analyzed demonstrating that the NIKS cells arose from the parental BC-1- Ep population. The odds of the NIKS cell line having the parental BC-1- Ep DNA fingerprint by random chance is  $4 \times 10^{-16}$ . The DNA fingerprints from three different sources of human

keratinocytes, ED-1-Ep, SCC4 and SCC13y are different from the BC-1-Ep pattern. This data also shows that keratinocytes isolated from other humans, ED-1-Ep, SCC4, and SCC13y, are unrelated to the BC-1-Ep cells or each other. The NIKS DNA fingerprint data provides an unequivocal way to identify the NIKS cell line.

**[0032]** Loss of p53 function is associated with an enhanced proliferative potential and increased frequency of immortality in cultured cells. The sequence of p53 in the NIKS cells is identical to published p53 sequences (GenBank accession number: M14695). In humans, p53 exists in two predominant polymorphic forms distinguished by the amino acid at codon 72. Both alleles of p53 in the NIKS cells are wild-type and have the sequence CGC at codon 72, which codes for an arginine. The other common form of p53 has a proline at this position. The entire sequence of p53 in the NIKS cells is identical to the BC-1-Ep progenitor cells. Rb was also found to be wild-type in NIKS cells.

**[0033]** Anchorage-independent growth is highly correlated to tumorigenicity *in vivo*. For this reason, the anchorage-independent growth characteristics of NIKS cells in agar or methylcellulose-containing medium were investigated. NIKS cells remained as single cells after 4 weeks in either agar- or methylcellulose-containing medium. The assays were continued for a total of 8 weeks to detect slow growing variants of the NIKS cells. None were observed.

**[0034]** To determine the tumorigenicity of the parental BC-1-Ep keratinocytes and the immortal NIKS keratinocyte cell line, cells were injected into the flanks of athymic nude mice. The human squamous cell carcinoma cell line, SCC4, was used as a positive control for tumor production in these animals. The injection of samples was designed such that animals received SCC4 cells in one flank and either the parental BC-1-Ep keratinocytes or the NIKS cells in the opposite flank. This injection strategy eliminated animal to animal variation in tumor production and confirmed that the mice would support vigorous growth of tumorigenic cells. Neither the parental BC-1-Ep keratinocytes (passage 6) nor the NIKS keratinocytes (passage 35) produced tumors in athymic nude mice.

**[0035]** NIKS cells were analyzed for the ability to undergo differentiation in both submerged culture and organotypic culture. Techniques for organotypic culture are described in detail in the examples. In particularly preferred embodiments, the organotypically cultured skin equivalents of the present invention comprise a dermal equivalent formed from collagen or a similar material and fibroblasts. The keratinocytes, for example NIKS cells or a combination of NIKS cells and cells from a patient are seeded onto the dermal equivalent and form an epidermal layer characterized by squamous differentiation following the organotypic culture process.

**[0036]** For cells in submerged culture, the formation cornified envelopes was monitored as a marker of squamous differentiation. In cultured human keratinocytes, early stages of cornified envelope assembly result in the formation of an immature structure composed of involucrin, cystatin-a and other proteins, which represent the innermost third of the mature cornified envelope. Less than 2% of the keratinocytes from the adherent BC-1-Ep cells or the NIKS cell

line produce cornified envelopes. This finding is consistent with previous studies demonstrating that actively growing, subconfluent keratinocytes produce less than 5% cornified envelopes. To determine whether the NIKS cell line is capable of producing cornified envelopes when induced to differentiate, the cells were removed from adherent culture and suspended for 24 hours in medium made semi-solid with methylcellulose. Many aspects of terminal differentiation, including differential expression of keratins and cornified envelope formation can be triggered *in vitro* by loss of keratinocyte cell-cell and cell-substratum adhesion. The NIKS keratinocytes produced as many as and usually more cornified envelopes than the parental keratinocytes. These findings demonstrate that the NIKS keratinocytes are not defective in their ability to initiate the formation of this cell type-specific differentiation structure.

**[0037]** To confirm that the NIKS keratinocytes can undergo squamous differentiation, the cells were cultivated in organotypic culture. Keratinocyte cultures grown on plastic substrata and submerged in medium replicate but exhibit limited differentiation. Specifically, human keratinocytes become confluent and undergo limited stratification producing a sheet consisting of 3 or more layers of keratinocytes. By light and electron microscopy there are striking differences between the architecture of the multilayered sheets formed in submerged culture and intact human skin. In contrast, organotypic culturing techniques allow for keratinocyte growth and differentiation under *in vivo*-like conditions. Specifically, the cells adhere to a physiological substratum consisting of dermal fibroblasts embedded within a fibrillar collagen base. The organotypic culture is maintained at the air-medium interface. In this way, cells in the upper sheets are air-exposed while the proliferating basal cells remain closest to the gradient of nutrients provided by diffusion through the collagen gel. Under these conditions, correct tissue architecture is formed. Several characteristics of a normal differentiating epidermis are evident. In both the parental cells and the NIKS cell line a single layer of cuboidal basal cells rests at the junction of the epidermis and the dermal equivalent. The rounded morphology and high nuclear to cytoplasmic ratio is indicative of an actively dividing population of keratinocytes. In normal human epidermis, as the basal cells divide they give rise to daughter cells that migrate upwards into the differentiating layers of the tissue. The daughter cells increase in size and become flattened and squamous. Eventually these cells enucleate and form cornified, keratinized structures. This normal differentiation process is evident in the upper layers of both the parental cells and the NIKS cells. The appearance of flattened squamous cells is evident in the upper epidermal layers and demonstrates that stratification has occurred in the organotypic cultures. In the uppermost part of the organotypic cultures the enucleated squames peel off the top of the culture. To date, no histological differences in differentiation at the light microscope level between the parental keratinocytes and the NIKS keratinocyte cell line grown in organotypic culture have been observed.

**[0038]** To observe more detailed characteristics of the parental (passage 5) and NIKS (passage 38) organotypic cultures and to confirm the histological observations, samples were analyzed using electron microscopy. Parental cells and the immortalized NIKS human keratinocyte cell line were harvested after 15 days in organotypic culture and sectioned perpendicular to the basal layer to show the extent of stratification. Both the parental cells and the NIKS cell line undergo extensive stratification in organotypic culture and form structures

that are characteristic of normal human epidermis. Abundant desmosomes are formed in organotypic cultures of parental cells and the NIKS cell line. The formation of a basal lamina and associated hemidesmosomes in the basal keratinocyte layers of both the parental cells and the cell line was also noted.

**[0039]** Hemidesmosomes are specialized structures that increase adhesion of the keratinocytes to the basal lamina and help maintain the integrity and strength of the tissue. The presence of these structures was especially evident in areas where the parental cells or the NIKS cells had attached directly to the porous support. These findings are consistent with earlier ultrastructural findings using human foreskin keratinocytes cultured on a fibroblast-containing porous support. Analysis at both the light and electron microscopic levels demonstrate that the NIKS cell line in organotypic culture can stratify, differentiate, and form structures such as desmosomes, basal lamina, and hemidesmosomes found in normal human epidermis.

## **B) Cryopreservation**

**[0040]** The present invention provides cryopreserved viable skin equivalents as defined in the claims. The cryopreserved skin equivalents are preferable storable at approximately -50 °C, -60 °C, -70 °C, -80 °C or colder for an extended period of time such as greater than 1, 2, 3, 4, 5 or 6 months and up to 12 or 24 months without a substantial loss of viability.

**[0041]** Preferably, all steps of the cryopreservation process prior to product packaging are performed aseptically inside a Class 100 biosafety cabinet in a Class 10,000 cleanroom. The cryopreservation process may comprise treating an organotypically cultured skin equivalent in a cryoprotectant solution. Certain examples of the present description are not limited to the use of any particular cryoprotectant. The cryoprotectant may be glycerol. The cryoprotectant may be provided in different concentrations in the cryoprotectant solution. The cryoprotectant may be provided in a solution comprising about 20% or 21% to about 70% of the solution by volume, and more preferably about 20% or 21% to about 45% of the solution by volume or 37.5% to 62.5% of the solution by volume, or most preferably from about 25% to 40% of the solution by volume or 42.5% to 57.5% of the solution by volume, depending on the temperature. The cryoprotectant solution may preferably comprise about 32.5% v/v or about 50% v/v cryoprotectant (e.g., glycerol). The cryoprotectant may be provided in a base medium solution. Suitable base medium solutions include, but are not limited to, DMEM, Ham's F-10, Ham's F-12, DMEM/F-12, Medium 199, MEM and RPMI. The base medium may form the remainder of the solution volume. The cryoprotectant solution may be buffered. Suitable buffers include, but are not limited to, HEPES, Tris, MOPS, and Trizma buffers. Buffering agents may be included at an amount to provide a buffered system in the range of pH 7.0 to 7.4. The cryoprotectant solution may be buffered with from about 5 mM to 15 mM HEPES, most preferably about 10 mM HEPES to a pH of about 7.0 to 7.4.

**[0042]** Treatment with the cryoprotectant solution may be conducted in a single step. By

"single step" it is meant that the cryoprotectant solution is not exchanged during the equilibration procedure as is common in the art. For example, the treatment step may be performed using a cryoprotectant solution with a defined concentration of cryoprotectant as opposed to a stepwise equilibration procedure where several media changes with increasing concentrations of cryoprotectant at each step. The treatment step may be conducted at a reduced temperature. The treatment step may be conducted at from about 2°C to 8°C, while in other examples, the treatment step may be conducted at room temperature, for example from about 15°C to 30°C. The skin equivalent may be incubated in the cryoprotectant solution for about 10 to 60 minutes, preferably from about 20 to 30 minutes.

**[0043]** The skin equivalent may be frozen following treatment with the cryoprotectant solution. Excess cryoprotectant solution may be removed from the skin equivalent prior to freezing, for example by aspirating the solution or moving the treated skin equivalent to a fresh vessel. Accordingly, the treated skin equivalent may be frozen by exposure to temperatures ranging from about -50 °C to -100 °C, and most preferably at about -80 °C. Treated skin equivalent may be simply placed in a bag or other vessel such as a culture dish and placed in a freezing unit such as a low temperature (e.g., -80 °C freezer) freezing unit. In contrast, it is common in the art to control the rate of freezing either by controlling the temperature in the freezing unit or by placing the tissue to be frozen in a container that allows control of the rate of decrease in temperature. The treated skin equivalent may be placed in a sterile culture vessel for freezing. The term "culture vessel" refers to any vessel of the type commonly used to culture cells or tissues and include circular, rectangular, and square dishes formed from a suitable material such as tissue culture plastic, polystyrene, polymers, plastics, glass, etc.

**[0044]** The cryopreserved skin equivalent may be packaged for long term storage. The skin equivalent may be placed in a bag or culture vessel as described above. The bag or culture vessel may be sealed, preferably heat sealed in a sterile bag (e.g., a plastic or polymer bag) to provide a primary package. The primary package is then sealed inside a secondary bag, for example a secondary plastic, foil, or Mylar bag. The cryopreserved tissues of the present description may preferably be stored at low temperature, from about -50 °C about -100 °C, preferably about -80 °C. The skin equivalents may be preferably stored from about 1, 2, 3, 4, 5 or 6 months and up to 12 or 24 months without a substantial loss of viability.

**[0045]** The present description provides a method of thawing a cryopreserved skin equivalent prior to application to a subject, comprising warming said cryopreserved skin equivalent and contacting said cryopreserved skin equivalent with an absorbent medium comprising a tissue compatible solution to allow removal of said cryoprotectant solution by diffusion. The cryopreserved skin equivalent in a suitable bag or culture vessel may be simply placed on a bench or table top and allowed to thaw. Thawing under controlled conditions as is common in the art is not necessary. Cryopreserved skin equivalent may be placed on an absorbent medium to remove thawed cryoprotectant solution from the skin equivalent. There is no limitation to the use a particular absorbent medium. Suitable absorbent media include, but are not limited to, Telfa pads, cellulosic pads (e.g., Whatman 1003-090 filter pads and Pall 70010 filter pads), gauze pads, and foam pads (e.g., Covidien 55544 hydrophilic foam pad). The

absorbent medium may be a Telfa pad. The absorbent medium may comprise a tissue-compatible solution. The tissue compatible solution may be a buffered solution. Suitable tissue compatible solutions include, but are not limited to, DMEM, Ham's F-10, Ham's F-12, DMEM/F-12, Medium 199, MEM and RPMI. Suitable buffers include, but are not limited to, HEPES, Tris, MOPS, and Trizma buffers. Buffering agents may be included at an amount to provide a buffered system in the range of pH 7.0 to 7.4.

**[0046]** The present description provides kits comprising a cryopreserved skin substitute, preferably provided in a package as described above. The kits may further comprise an absorbent medium, and a tissue compatible solution.

**[0047]** The present description provides a process for forming an organotypically cultured skin equivalent and freezing the skin equivalent in the same culture vessel. The culture vessel may comprise an insert movable between upper and lower positions in the culture dish. The insert preferably comprises a bottom surface which is a porous membrane. In use, the vessel is filled with the appropriate culture medium and a dermal equivalent is formed on the porous membrane. The dermal equivalent is then seeded with keratinocytes (e.g., NIKS cells) in the presence of the appropriate culture medium. At the appropriate time, an air interface is created by lowering the level of culture medium in the vessel and the culture is continued until the stratified skin equivalent is formed. The culture medium is then removed from the vessel and the insert is lowered to the lower position. The cryoprotectant solution is added for treatment and then removed, and the vessel is then transferred to the freezing unit.

### **C) Therapeutic Uses**

**[0048]** It is contemplated that the cryopreserved skin equivalents of the present invention as defined in the claims may be used therapeutically. The cryopreserved skin substitute may be used in wound closure and burn treatment applications. The use of autografts and allografts for the treatment of burns and wound closure is described in Myers et al., A. J. Surg. 170(1):75-83 (1995) and U.S. Pat. Nos. 5,693,332; 5,658,331; and 6,039,760. The skin equivalents may be used in conjunction with dermal replacements such as DERMAGRAFT or INTEGRA. Accordingly, the present description provides methods for wound closure, including ulcers or wounds caused by burns, comprising providing a skin equivalent and a patient suffering from a wound and treating the patient with the skin equivalent under conditions such that the wound is closed.

**[0049]** The skin equivalents may be utilized to treat chronic skin wounds. Chronic skin wounds (e.g., venous ulcers, diabetic ulcers, pressure ulcers) are a serious problem. The healing of such a wound often takes well over a year of treatment. Treatment options currently include dressings and debridement (use of chemicals or surgery to clear away necrotic tissue), and/or antibiotics in the case of infection. These treatment options take extended periods of time and high levels of patient compliance. As such, a therapy that can increase a practitioner's success in healing chronic wounds and accelerate the rate of wound healing would meet an unmet



need in the field. Accordingly, the present description contemplates treatment of skin wounds with cryopreserved skin equivalents. Skin equivalents may be topically applied to wounds. Cryopreserved skin equivalents may be used for application to partial thickness wounds. Cryopreserved skin equivalents may be used to treat full thickness wounds. Cryopreserved skin equivalents may be used to treat numerous types of internal wounds, including, but not limited to, internal wounds of the mucous membranes that line the gastrointestinal tract, ulcerative colitis, and inflammation of mucous membranes that may be caused by cancer therapies. Skin equivalents expressing host defense peptides or pro-angiogenic factors may be used as a temporary or permanent wound dressing.

**[0050]** The cells may be engineered to provide additional therapeutic agents to a subject. There is no limitation to the delivery of any particular therapeutic agent. Indeed, it is contemplated that a variety of therapeutic agents may be delivered to the subject, including, but not limited to, enzymes, peptides, peptide hormones, other proteins, ribosomal RNA, ribozymes, small interfering RNA (siRNA) micro RNA (miRNA), and antisense RNA. The agents may be host defense peptides such as human beta-defensin 1, 2, or 3 or cathelicidin or other proteins such as VEGF and HIF-1 $\alpha$ , see, e.g., U.S. Pat. Nos. 7,674,291; 7,807,148; 7,915,042; 7,988,959; and 8,092,531. These therapeutic agents may be delivered for a variety of purposes, including but not limited to the purpose of correcting genetic defects. The therapeutic agent may be delivered for the purpose of detoxifying a patient with an inherited inborn error of metabolism (e.g., aminoacidopathesis) in which the skin equivalent serves as wild-type tissue. It is contemplated that delivery of the therapeutic agent corrects the defect. The cells may be transfected with a DNA construct encoding a therapeutic agent (e.g., insulin, clotting factor IX, erythropoietin, etc.) and skin equivalents prepared from transfected cells are administered to the subject. The therapeutic agent is then delivered to the patient's bloodstream or other tissues from the graft. The nucleic acid encoding the therapeutic agent may be operably linked to a suitable promoter. There is no limitation to the use of any particular promoter. Indeed, the use of a variety of promoters is contemplated, including, but not limited to, inducible, constitutive, tissue-specific, and keratinocyte-specific promoters. The nucleic acid encoding the therapeutic agent may be introduced directly into the keratinocytes (i.e., by electroporation, calcium phosphate co-precipitation, or liposome transfection). The nucleic acid encoding the therapeutic agent may be provided as a vector and the vector is introduced into the keratinocytes by methods known in the art. The vector may be an episomal vector such as a replicating plasmid. The vector may integrate into the genome of the keratinocytes. Examples of integrating vectors include, but are not limited to, retroviral vectors, adeno-associated virus vectors, non-replicating plasmid vectors and transposon vectors.

## **EXPERIMENTAL**

**[0051]** The following examples are provided.

**[0052]** In the experimental disclosure which follows, the following abbreviations apply: eq (equivalents); M (Molar); mM (millimolar);  $\mu$ M (micromolar); N (Normal); mol (moles); mmol

(millimoles);  $\mu\text{mol}$  (micromoles); nmol (nanomoles); g (grams); mg (milligrams);  $\mu\text{g}$  (micrograms); ng (nanograms); 1 or L (liters); ml or mL (milliliters);  $\mu\text{l}$  or  $\mu\text{L}$  (microliters); cm (centimeters); mm (millimeters);  $\mu\text{m}$  (micrometers); nm (nanometers);  $^{\circ}\text{C}$  (degrees Centigrade); U (units), mU (milliunits); min. (minutes); sec. (seconds); % (percent); kb (kilobase); bp (base pair); PCR (polymerase chain reaction); BSA (bovine serum albumin); CFU (colony forming units); kGy (kiloGray); PVDF (polyvinylidene fluoride); BCA (bicinchoninic acid); SDS-PAGE (sodium dodecyl sulfate polyacrylamide gel electrophoresis).

#### **Example 1 (not according to the invention)**

**[0053]** StrataGraft® skin tissue is a living, full-thickness, allogeneic human skin substitute that reproduces many of the structural and biological properties of normal human skin. StrataGraft® skin tissue contains both a viable, fully-stratified epidermal layer derived from NIKS® cells, which are a consistent and well-characterized source of pathogen-free human keratinocyte progenitors, and a dermal layer containing normal human dermal fibroblasts (NHDF) embedded in a collagen-rich matrix. StrataGraft® skin tissue possesses excellent tensile strength and handling characteristics that enable it to be meshed, stapled, and sutured similarly to human skin grafts. StrataGraft® also exhibits barrier function comparable to that of intact human skin and is capable of delivering bioactive molecules for wound bed conditioning and tissue regeneration. The physical and biological characteristics of StrataGraft skin tissue make it ideal for the treatment of a variety of skin wounds.

**[0054]** The manufacturing process for StrataGraft® skin tissue encompasses three sequential cell and tissue culture processes. In Stage I of the manufacturing process, NIKS keratinocytes are expanded in monolayer cell culture. Concurrent with the NIKS keratinocyte culture in Stage I, NHDF are expanded in monolayer culture and combined with purified type I collagen and culture medium and allowed to gel to form the cellularized dermal equivalent (DE). Alternatively, NHDF are seeded into Transwell inserts and allowed to proliferate and secrete and assemble extracellular matrix molecules into a simplified dermal equivalent. In Stage II, NIKS keratinocytes are seeded onto the surface of the DE and cultured under submerged conditions for two days to promote complete epithelialization of the DE surface. The tissue is then lifted to the air-liquid interface in Stage III, where it is maintained for 18 days in a controlled, low humidity environment to promote tissue maturation. The skin equivalents are generally prepared as described in U.S. Pat. Nos. 7,674,291; 7,807,148; 7,915,042; 7,988,959; and 8,092,531.

#### **Example 2**

**[0055]** This example describes improved cryopreservation methods for human skin equivalents. The production process is unchanged from the current method described above. All tissues in the lot are fed with fresh medium and incubated overnight prior to

cryopreservation. Prior to cryopreservation, media samples from all tissues in each lot are tested for sterility. The remaining tissues in each lot are cryopreserved as follows.

Parameter	Operating Range
Cryoprotectant formulation	50% (v/v) glycerol
	DMEM (IX)
	10 mM HEPES (pH 7.0 to 7.4)
Pre-freeze cryoprotectant incubation temperature	2-8 °C
Pre-freeze cryoprotectant incubation time	20-30 minutes
Freeze method	Direct transfer to -80 °C freezer
Storage temperature	-70 to -90 °C
Shipping conditions	Overnight delivery on dry ice

### Cryopreservation Process Description

**[0056]** All steps of the cryopreservation process prior to the final product packaging step are performed aseptically inside a Class 100 biosafety cabinet in a Class 10,000 cleanroom.

Step 1- Pre-cool 100 mm culture dishes containing 20 ml of cryoprotectant solution to 2-8°C on a stainless steel cold treatment surface inside biosafety cabinet. Temperature of cold treatment surface is maintained at 2-8°C for several hours by contact with frozen gel packs submerged in water.

Step 2- Transfer Transwells containing StrataGraft® tissues into individual dishes containing pre-cooled cryoprotectant solution. Incubate tissues 20-30 minutes in cryoprotectant on the cold treatment surface.

Step 3- Transfer treated StrataGraft® tissues to new sterile 100 mm culture dishes containing final product label so that the tissue rests on the bottom of the culture dish and return tissues back to the cold treatment surface. Excess cryoprotectant is allowed to drain from the skin equivalent to provide a treated skin equivalent that is substantially free of excess cryoprotectant on the exterior surfaces of the skin equivalent.

Step 4- Heat-seal 100 mm culture dishes in clear, sterile bags. Place primary package into secondary Mylar bag, heat-seal, and transfer packaged tissues to cold storage container until all tissues are packaged.

Step 5- Remove cold storage container with packaged StrataGraft® tissues from cleanroom and transfer tissues to an ultralow freezer (-75°C to -80°C). Place tissues in a pre-cooled rack in the freezer that allows unrestricted airflow to the top and bottom of the packaged tissues to ensure uniform and rapid cooling. Leave tissues undisturbed overnight during the freezing process.

**[0057]** Tissues are placed into quarantine storage at -70 to -90°C pending results of lot release testing. A representative tissue from each lot of cryopreserved StrataGraft® skin tissue is tested using a panel of Quality Control SOP that have historically been used for lot release testing of StrataGraft® tissue.

**[0058]** Stratatech has established and qualified a panel of lot release assays that are used to characterize StrataGraft® skin tissue. A subset of these lot release assays has been used to monitor and evaluate the impact that changes to the storage conditions may have on key biological and structural characteristics of StrataGraft® skin tissue (e.g., barrier function, viability, and histological appearance). Although transient minor changes in the histological appearance of StrataGraft® tissue are generally observed following cryopreservation, the histological architecture normalizes after being reintroduced into organotypic culture for several days, indicating that the viable cells in the basal layer of StrataGraft® are able to proliferate and reproduce the epidermal layer after cryopreservation. The systematic evaluation of cryoprotectant concentration, incubation time, incubation temperature, freezing rate, and storage conditions has enabled Stratatech to identify a cryopreservation process that enables long-term storage of StrataGraft® skin tissue with consistent and defined quality and that meets the specifications that define StrataGraft® skin tissue.

**[0059]** The following parameters were systematically evaluated during development of the StrataGraft cryopreservation process.

- Cryoprotectant composition
- Pre-freeze cryoprotectant incubation temperature
- Pre-freeze cryoprotectant incubation time
- Number of steps required during cryoprotectant incubation
- Freezing rate
- Final product packaging
- Storage temperature
- Shipping conditions
- Thaw temperature and time
- Post-thaw cryoprotectant diffusion mediator
- Post-thaw incubation solution
- Post-thaw incubation temperature
- Post-thaw incubation time

**[0060]** As anticipated, many of these individual parameters interact to influence the properties of cryopreserved tissue. For example, it is not possible to optimize cryoprotectant concentration without also taking into account the cryoprotectant incubation time and temperature. Likewise, post-thaw incubation temperature influences the allowable range of

post-thaw incubation times. During development of the cryopreservation process, a range of acceptable values for each of the individual parameters was identified and used to define the final combination of cryoprotectant formulation, pre-freeze incubation time, pre-freeze incubation temperature, freeze rate, and storage condition. The operating parameters of the cryopreservation process as developed for cryopreservation and storage of StrataGraft® skin tissue are listed above.

**[0061]** Glycerol (glycerin) was identified as the most desirable cryoprotectant for StrataGraft® tissue. The glycerol used in the cryoprotectant formulation is synthetic, USP-grade material that undergoes additional testing for endotoxin prior to release for use. In addition to glycerol, the cryoprotectant solution contains Dulbecco's Modified Eagle Medium (DMEM) and 10 mM HEPES to maintain pH of 7.0 to 7.4 at ambient atmospheric conditions. DMEM was chosen as the base for the cryopreservation solution because it is already a component of the culture medium used to prepare StrataGraft® skin tissue. HEPES is a well-characterized buffering agent that maintains the pH of the cryopreservation solution outside of a CO<sub>2</sub> environment.

**[0062]** A series of studies were performed to determine an appropriate glycerol concentration for use in the cryoprotectant solution. Glycerol concentrations tested ranged from 16.25% to 65%. In some cases (which are not according to the present invention), the concentration of glycerol was gradually increased in two or three steps by incubation in a series of solutions with increasing glycerol concentration. Initial studies to demonstrate the feasibility of cryopreserving StrataGraft® tissue used a three-step process in which the glycerol concentration was sequentially increased to 16.25%, then 32.5%, and finally to 65% by incubation for 15-20 minutes in each glycerol solution, while gradually reducing the temperature at each incubation step (16.25% glycerol incubation at room temperature, 32.5% at 2-8°C, and 65% at -20°C). Finally, tissues were frozen to -140°C at -15 °C/min in a controlled-rate freezer. Cryopreserved tissues were transferred to an ultra-cold freezer (-80°C), where they were maintained until thawed for analysis. While this three-step process could be used to preserve the viability and histological architecture of StrataGraft® tissues during cryopreserved storage, the complexity of having multiple steps with different solutions performed at different temperatures introduces the opportunity for error and is not amenable to process scale-up that would be required for commercialization.

**[0063]** Analysis of tissues cryopreserved with this three-step method revealed that, after thawing and dilution of the cryoprotectant, the tissues exhibited viability (assessed by the ability of viable cells in the tissue to convert MTT to its formazan product), histological architecture, and barrier function comparable to tissues that were not cryopreserved. The tissues maintained high levels of viability following re-introduction into organotypic culture for up to nine days after thawing, demonstrating that the metabolic activity detected shortly after thawing was not just residual enzymatic activity.

**[0064]** During these initial studies, it was found that the concentration of glycerol in cryopreserved tissues could be reduced after thawing either by incubating the tissues in a series of solutions with decreasing glycerol concentration, or by placing the tissues in a media

reservoir with a filter pad just below the tissue to moderate the diffusion of glycerol. Based on these results, subsequent studies primarily used the single-step approach of incubating thawed tissues in culture medium with a pad to moderate glycerol diffusion. See Figure 7.

**[0065]** After demonstrating the feasibility of cryopreserving StrataGraft® skin tissue using the three-step method, we used this method as a benchmark against which to compare simplifications to the cryopreservation method. These studies examined reducing the number of steps required to reach the same final glycerol concentration (65%) or evaluating final glycerol concentrations that were lower than that used in the initial studies. Based on these studies, it was determined that glycerol concentrations as low as 32.5% could be used to reproducibly maintain the viability and histological architecture of StrataGraft® tissue during cryopreserved storage. In contrast, a final glycerol concentration of 16.25% in the cryoprotectant solution did not support maintenance of viability in frozen tissues. See Figure 5. By evaluating a range of glycerol concentrations, it was determined that a cryoprotectant solution containing 50% glycerol reproducibly supported cryopreservation of StrataGraft® skin tissue and provided a margin of error above somewhat lower glycerol concentrations (e.g., 32.5%) that also supported efficient cryopreservation.

***Pre-freeze cryoprotectant incubation temperature, incubation time, and number of incubation steps***

**[0066]** In addition to the concentration of glycerol in the cryoprotectant solution, three other factors affecting treatment of StrataGraft® tissue with cryoprotectant prior to cryopreservation are: 1) the pre-freeze cryoprotectant incubation time, 2) the temperature at which the tissues are treated with cryoprotectant, and 3) the number of steps required to reach the final glycerol concentration. As described above, initial feasibility studies with the three-step process involved sequentially incubating the tissues in solutions with increasing glycerol concentrations at successively lower temperatures for 15-20 minutes at each step. As stated above, a simpler cryopreservation process is preferred to avoid the complexity of having multiple steps with different solutions performed at different temperatures to reduce the opportunity for error and facilitate process scale-up. Toward this goal, the need for stepwise increase in cryoprotectant concentration and stepwise reduction in temperature during the pre-freeze cryoprotectant equilibration phase was evaluated.

**[0067]** In a series of studies performed in conjunction with evaluation of different cryoprotectant concentrations, it was determined that StrataGraft® tissues treated with cryoprotectant solutions containing 32.5, 50, or 65% glycerol in a single step at 2-8°C for as little as 15 minutes and as long as 60 minutes were all able to withstand cryopreservation with minimal loss of viability or epidermal architecture. Although no decline in tissue performance was observed with cryoprotectant treatment times up to 60 minutes, relatively short glycerol treatment times (20 to 30 minutes) were chosen in order to minimize any potential adverse effects of prolonged exposure to cryoprotectant prior to freezing.

### ***Freezing rate***

**[0068]** As described above, the initial feasibility studies utilized a controlled-rate freezer to freeze tissues at a rate of  $-15^{\circ}\text{C}/\text{min}$  after equilibration with cryoprotectant. However, the use of a controlled-rate freezer would impose significant additional costs and is not amenable to process scale-up. Historically, cryopreservation of human cells has been accomplished at a more moderate rate of approximately  $-1^{\circ}\text{C}/\text{min}$  without the use of controlled-rate freezers. This is routinely accomplished by placing vials of cells in an ultra-cold freezer in an insulated box or container designed to moderate the cooling rate to approximately  $-1^{\circ}\text{C}/\text{min}$ .

**[0069]** Process development studies were designed to determine whether StrataGraft® tissues could be cryopreserved without using a controlled-rate freezer. These studies demonstrated that tissues frozen by direct transfer to an ultra-cold freezer (approximately  $-80^{\circ}\text{C}$ ) after treatment with cryoprotectant performed as well as tissues frozen in a controlled-rate freezer. See Figure 1.

**[0070]** Temperature monitoring studies were performed to track the temperature of tissues during the cryopreservation process. Following packaging in an inner sterile bag and an outer Mylar bag, StrataGraft® tissues are transferred to a pre-cooled rack inside of an ultra-cold freezer. Each tissue is placed in a separate slot in the freezer rack, with ample room above and below the tissue to allow unrestricted airflow during the freezing process. Using temperature monitoring probes positioned within culture dishes packaged as described above and loaded into freezer racks in this configuration, the temperature rapidly decreases to approximately  $-50^{\circ}\text{C}$  within the first 15 minutes, further cools to approximately  $-65^{\circ}\text{C}$  by 30 minutes and reaches a final temperature of approximately  $-80^{\circ}\text{C}$  after three hours. There is no significant difference in the temperature profiles between tissues placed in the top, middle, and bottom positions of the freezer rack. See Figure 2.

### ***Final product packaging***

**[0071]** In initial studies, tissues were frozen in contact with a layer of cryoprotectant solution after incubation with cryoprotectant. Although tissues cryopreserved in this manner exhibited good post-thaw properties, rapid thawing of tissues frozen in contact with this layer of cryoprotectant required incubation for several minutes in a  $35\text{-}39^{\circ}\text{C}$  water bath, which would be difficult to implement and standardize in a surgical suite. It was subsequently determined that contact with the cryoprotectant solution was not required after tissues had been treated with cryoprotectant. This enabled development of the final product configuration in which tissues are transferred to an empty sterile 100 mm culture dish after treatment with cryoprotectant, where they are frozen in contact with the bottom of the empty dish rather than being frozen in contact with a layer of cryoprotectant solution.

**[0072]** To maintain the sterility of cryopreserved tissues, the 100 mm culture dishes containing cryoprotectant-treated tissues are aseptically packaged and heat-sealed inside of a sterile polyethylene sample bag. The inner bag is then heat-sealed inside a puncture-resistant, food grade, metalized polyester/polyethylene bag, which protects the packaged tissues from light, moisture, and provides a barrier to CO<sub>2</sub> vapor during shipment on dry ice. The stability and comparability studies described below utilized tissues packaged and cryopreserved in this configuration.

### ***Storage temperature***

**[0073]** Cryopreservation of viable skin equivalents enables burn centers to have ready access to this cell-based regenerative medicine therapeutic for burns and other indications that require rapid intervention. Optimally, major burn centers would be able to maintain an inventory of the product for use without the need to schedule a delivery on a case by case basis. In early feasibility studies, cryopreserved tissues were stored at -196°C in a vapor-phase nitrogen freezer. Since burn centers do not typically have liquid nitrogen storage capabilities, cryopreservation procedures were developed that permit storage of tissue for at least six months in ultra-cold freezers (-60°C to -90°C), which are readily available in blood and tissue banks at most hospitals and trauma centers. Results of these experiments, demonstrated that while tissues stored at -50°C exhibited significant losses of viability over the course of several weeks, tissues stored at -80°C retained levels of viability comparable to tissues that had been stored in nitrogen vapor. See Figure 3. These results were obtained with several independent lots of cryopreserved tissue, confirming the reproducibility of this finding. As described in the *Stability of Cryopreserved Tissues* section below, analysis of tissues stored at -80°C for six months revealed no significant loss of viability or changes to epidermal architecture during storage.

### ***Shipping conditions***

**[0074]** Cryopreserved StrataGraft® skin tissue will be shipped to clinical sites on dry ice for next morning delivery via commercial courier such as FedEx or UPS. The shipping container (Freezetherm FT29, Laminar Medica) is a validated dry ice shipping box that holds sufficient dry ice to maintain the cryopreserved tissues at <-75°C for at least 72 hours at ambient temperatures of up to 35°C to account for possible delays in delivery. Experimental data indicates that storage of cryopreserved tissues in the dry ice shipping container for >48 hours does not have any detectable adverse effect on tissue viability or histological architecture. See Figure 4. Following receipt at the clinical site, cryopreserved StrataGraft tissues will be stored in an ultra-cold freezer (e.g., -60°C to -90°C) until use.

### ***Pre-operative preparation of cryopreserved tissue***



**[0075]** Prior to clinical use, cryopreserved StrataGraft tissue is thawed and incubated briefly on pads saturated with culture medium to remove residual cryoprotectant. Due to the geometry of the tissue, the thaw phase is rapid. After tissues are thawed, Transwell inserts containing the tissue are aseptically transferred to the sterile field and placed in sterile dishes containing absorbent pads saturated with culture medium. As described below, the timing of the post-thaw incubation phase is flexible enough to accommodate delays that could be reasonably anticipated during clinical use.

#### ***Thaw temperature and time***

**[0076]** As described in the *Final product packaging* section above, tissues are cryopreserved in a culture dish without a layer of cryoprotectant solution, which allows the tissues to be thawed rapidly at ambient temperature simply by placing the package onto a bench or table. Precise control over the thaw temperature and time is not required, as experimental data shows that tissues thawed for varying times at temperatures ranging from 22°C to 40°C exhibit similar post-thaw properties.

#### ***Post-thaw incubation solution***

**[0077]** Buffered post-thaw incubation solutions work better than unbuffered solutions. Tissues incubated in simple unbuffered salt solutions (lactated Ringer's or normal saline) do not survive as well as tissues incubated in culture media-based solutions. Stratatech's SM01 culture media (StrataLife series) or commercially available DMEM/F12 media buffered with HEPES are preferred. See Figure 6.

#### ***Post-thaw incubation temperature***

**[0078]** Initial development studies found that buffered culture media were able to support tissue viability following post-thaw incubation at both 37 °C and room temperature. However, warm post-thaw incubation temperatures (37 °C) work better than cooler temperatures (20-25 °C) for sub-optimal incubation solutions. See Figure 6. Starting the post-thaw incubation on a pad containing 37 °C media that slowly cools to room temp over the course of 15-30 minutes also works well. Higher temperature seems to be most important in the first few minutes after thaw.

**[0079]** Later development studies demonstrated that buffered media solutions pre-warmed to only room temperature were comparable to media warmed to 37 °C in their ability to support the properties of thawed tissues. This was true for tissues held on either Telfa pads or Whatman pads. See Figure 8

**Post-thaw incubation time**

[0080] Tissues can be left on the media saturated pad for times ranging from 15 min to 4 hr at 20-25°C or up to 2 hr at 40°C with no significant effect on tissue viability.

**STABILITY STUDY**

**Stability of cryopreserved tissues**

[0081] Although many of the studies described above analyzed tissues that had been stored in a frozen state for only a few days or weeks, it is widely accepted that the majority of damage to cryopreserved cells and tissues occurs during the freezing and thawing stages, with relatively little loss of cellular viability taking place during long-term storage at reduced temperatures. Long-term storage results show that cryopreserved StrataGraft® skin tissue maintains high levels of viability and histological architecture after at least 12 months of storage at ultracold temperatures. Analysis of tissues produced and cryopreserved using the cryopreservation process described above indicates that tissues cryopreserved with this process maintain key biological, structural, and physical properties during storage for at least 12 months at ultracold temperatures. See Figure 9.

**Example 3**

[0082] This example describes improved cryopreservation methods for human skin equivalents utilizing a pre-freeze treatment step with cryopreservation solutions containing 32.5% or 50% glycerol at room temperature. The production process is unchanged from the current method described previously. At the end of the production process, the tissues are treated and cryopreserved as follows.

Parameter	Operating Range
Cryoprotectant formulation	32.5% (v/v) glycerol
	DMEM (IX)
	10 mM HEPES (pH 7.0 to 7.4)
	or
	50% (v/v) glycerol
	DMEM (IX)
	10 mM HEPES (pH 7.0 to 7.4)
Pre-freeze cryoprotectant incubation temperature	Room temperature

Parameter	Operating Range
Pre-freeze cryoprotectant incubation time	15-45 minutes
Freeze method	Direct transfer to -80 °C freezer
Storage temperature	-70 to -90 °C
Shipping conditions	Overnight delivery on dry ice

### Cryopreservation Process Description

**[0083]** All steps of the cryopreservation process prior to the final product packaging step are performed aseptically inside a Class 100 biosafety cabinet in a Class 10,000 cleanroom.

Step 1- Dispense 20 ml of cryoprotectant solution to 100 mm culture dishes.

Step 2- Transfer Transwells containing StrataGraft® tissues into individual dishes containing cryoprotectant solution. Incubate tissues 15-45 minutes in cryoprotectant.

Step 3- Transfer treated StrataGraft® tissues to new sterile 100 mm culture dishes containing final product label so that the tissue rests on the bottom of the culture dish. Excess cryoprotectant is allowed to drain from the skin equivalent to provide a treated skin equivalent that is substantially free of excess cryoprotectant on the exterior surfaces of the skin equivalent.

Step 4- Heat-seal 100 mm culture dishes in clear, sterile bags. Place primary package into secondary Mylar bag and heat-seal.

Step 5- Remove the packaged StrataGraft® tissues from cleanroom and transfer tissues to an ultralow freezer (-75°C to -80°C). Place tissues in a pre-cooled rack in the freezer that allows unrestricted airflow to the top and bottom of the packaged tissues to ensure uniform and rapid cooling. Leave tissues undisturbed overnight during the freezing process.

**[0084]** Cryopreserved tissues were thawed at room temperature for 10 minutes, transferred to a hold chamber containing Telfa pads saturated with 40 ml of HEPES-buffered culture medium that had been warmed to room temperature, and held at RT for 15 to 20 minutes. Tissues were transferred to a culture dish containing 90 ml of SM01 medium and returned to culture overnight. Tissues were analyzed for viability after overnight re-culture. Tissues pre-treated with 32.5% glycerol at room temperature for 15 to 45 minutes had acceptable post-thaw viability. Tissues treated with 50% glycerol at room temperature for 15 minutes also had acceptable viability; however, tissues treated with 50% glycerol at room temperature for 45 minutes had unacceptable viability. See Figure 10.

### Example 4

**[0085]** MTT assays are preferably conducted as follows. Samples are excised from the skin tissue using an 8 mm biopsy punch. The samples are transferred to 0.3 ml MTT Assay Medium (1 mg/ml MTT reagent in Ham's F-12) in a 24-well plate that has been pre-warmed to 37 °C / 5% CO<sub>2</sub>. The samples are incubated for 85-95 minutes at 37 °C / 5% CO<sub>2</sub>. The samples are blotted and transferred to 2 ml isopropanol. The samples are thoroughly mixed to completely extract the purple formazan product. 200 µl in triplicate of each extract is aliquoted into a 96-well plate, using isopropanol as a blank. The absorbance (550 nm) is measured in a spectrophotometer.

tissues to ensure uniform and rapid cooling. Leave tissues undisturbed overnight during the freezing process.

**[0086]** Cryopreserved tissues were thawed at room temperature for 10 minutes, transferred to a hold chamber containing Telfa pads saturated with 40 ml of HEPES-buffered culture medium that had been warmed to room temperature, and held at RT for 15 to 20 minutes. Tissues were transferred to a culture dish containing 90 ml of SM01 medium and returned to culture overnight. Tissues were analyzed for viability after overnight re-culture. Tissues pre-treated with 32.5% glycerol at room temperature for 15 to 45 minutes had acceptable post-thaw viability. Tissues treated with 50% glycerol at room temperature for 15 minutes also had acceptable viability; however, tissues treated with 50% glycerol at room temperature for 45 minutes had unacceptable viability. See Figure 10.

#### **Example 4**

**[0087]** MTT assays are preferably conducted as follows. Samples are excised from the skin tissue using an 8 mm biopsy punch. The samples are transferred to 0.3 ml MTT Assay Medium (1 mg/ml MTT reagent in Ham's F-12) in a 24-well plate that has been pre-warmed to 37 °C / 5% CO<sub>2</sub>. The samples are incubated for 85-95 minutes at 37 °C / 5% CO<sub>2</sub>. The samples are blotted and transferred to 2 ml isopropanol. The samples are thoroughly mixed to completely extract the purple formazan product. 200 µl in triplicate of each extract is aliquoted into a 96-well plate, using isopropanol as a blank. The absorbance (550 nm) is measured in a spectrophotometer.

## **REFERENCES CITED IN THE DESCRIPTION**

This list of references cited by the applicant is for the reader's convenience only. It does not form part of the European patent document. Even though great care has been taken in compiling the references, errors or omissions cannot be excluded and the EPO disclaims all liability in this regard.

**Patent documents cited in the description**

- [US7674291B](#) [0026] [0050] [0054]
- [US7807148B](#) [0026] [0050] [0054]
- [US7915042B](#) [0026] [0050] [0054]
- [US7988959B](#) [0026] [0050] [0054]
- [US8092531B](#) [0026] [0050] [0054]
- [US5968546A](#) [0028]
- [US5693332A](#) [0028] [0048]
- [US4485096A](#) [0028]
- [US6039760A](#) [0028] [0048]
- [US5536656A](#) [0028]
- [US5989837A](#) [0028]
- [US5658331A](#) [0048]

**Non-patent literature cited in the description**

- **DE BREE** Genomics-based Drug Data Report and Regenerative Therapy, 2001, vol. 2, 177-96 [0002]
- BioMedGPS-SmarTRAK market analysis, 2012, [0002]
- **KAGAN et al.** Clin Lab Med, 2005, vol. 25, 587-605 [0025]
- **AUGER et al.** In Vitro Cell. Dev. Biol. - Animal, vol. 36, 96-103 [0028]
- **ASBILL et al.** Pharm. Research, 2000, vol. 17, 91092-97 [0028]
- **MEANA et al.** Burns, 1998, vol. 24, 621-30 [0028]
- **BADEN** In Vitro Cell. Dev. Biol., 1987, vol. 23, 3205-213 [0028]
- **BOUCAMP et al.** J. cell. Boil., 1988, vol. 106, 761-771 [0028]
- **ALLEN-HOFFMANN et al.** J. Invest. Dermatol., 2000, vol. 114, 3444-455 [0028]
- **MYERS et al.** A. J. Surg., 1995, vol. 170, 175-83 [0048]

## **Patentkrav**

1. Fremgangsmåde til kryokonservering af en organotypisk dyrket hudækvivalent til opretholdelse af levedygtigt væv, omfattende:
  - a) at behandle en organotypisk dyrket hudækvivalent i en kryobeskyttende opløsning i et enkelt trin, hvor den organotypisk dyrkede hudækvivalent omfatter lagdelt pladeepitel på et dermalt lag omfattende fibroblast, og hvor det kryobeskyttende middel er tilvejebragt i en opløsning omfattende 21 til 70 volumen%, fortrinsvis 21 til 45 volumen% eller 37,5 til 62,5 volumen%, mere foretrukket 25 til 40 volumen% eller 42,5 til 57,5 volumen% af opløsningen, og det kryobeskyttende middel er glycerol;
  - b) at separere den behandlede organotypisk dyrkede hudækvivalent fra overskydende kryobeskyttende opløsning og emballere den behandlede organotypisk dyrkede hudækvivalent i fravær af yderligere kryobeskyttende middel til tilvejebringelse af en emballeret hudækvivalent, hvor emballeringen fortrinsvis endvidere omfatter at indeslutte den kryokonserverede hudækvivalent i en steril pose og indeslutte den sterile pose i en anden pose; og
  - c) at fryse den emballerede organotypisk dyrkede hudækvivalent til tilvejebringelse af en kryokonserveret hudækvivalent.
2. Fremgangsmåde ifølge krav 1, hvor behandlingen af en organotypisk dyrket hudækvivalent i en kryobeskyttende opløsning i et enkelt trin er ved en temperatur på 2° C til 8° C eller 15° C til 30° C.
3. Fremgangsmåde ifølge krav 1, hvor frysningen endvidere omfatter frysning ved -80° C eller direkte udsættelse for temperaturer i området fra -50° C til -100° C.
4. Fremgangsmåde ifølge krav 1, hvor de organotypisk dyrkede hudækvivalenter omfatter NIKS-celler (nær-diploide immortaliserede keratinocytter), hvor NIKS-cellerne fortrinsvis omfatter en eksogen nukleinsyresekvens, der koder for et eksogent polypeptid.
5. Fremgangsmåde ifølge krav 1, hvor hudækvivalenten bibeholder levedygtighed efter optøning, hvor den kryokonserverede hudækvivalent fortrinsvis har en A<sub>550</sub> på mindst 50% af en referencehudækvivalent som bestemt ved

hjælp af et MTT-assay.

- 5 **6.** Fremgangsmåde til optøning af en kryokonserveret hudækvivalent, der er fremstillet ved hjælp af fremgangsmåden ifølge krav 1, til påføring på en patient, omfattende:
- a) at opvarme den kryokonserverede hudækvivalent, hvor opvarmningen fortrinsvis omfatter udsættelse for rumtemperatur på anvendelsesstedet; og
- b) at bringe den kryokonserverede hudækvivalent i kontakt med en diffusionsmediator omfattende en vævskompatibel opløsning, fortrinsvis en bufret vævskompatibel opløsning, til fjernelse af den kryobeskyttende opløsning ved diffusion, og hvor diffusionsmediatoren er valgt fra gruppen bestående af et absorberende medium, hvor det absorberende medium fortrinsvis er valgt fra gruppen bestående af
- 10 - Telfa-pads, skumpads, gaze pads og cellulose pads indeholdende det vævskompatible medium,
- 15 - en membran og
- en dialysepose.
- 20 **7.** Kryokonserveret hudækvivalent fremstillet ved hjælp af fremgangsmåden ifølge krav 1.
- 8.** System omfattende hudækvivalenten ifølge krav 7, der er anbragt på et absorberende medium.
- 25 **9.** Emballeret kryokonserveret hudækvivalent ifølge krav 7 til anvendelse ved behandling af sår, hvor hudækvivalenten er til påføring på et sår under betingelser, hvor hudækvivalenten kommer i kontakt med såret.
- 30 **10.** Kit omfattende en kryokonserveret hudækvivalent ifølge krav 7, en diffusionsmediator og en vævskompatibel opløsning.
- 11.** Kit ifølge krav 10, hvor den kryokonserverede hudækvivalent er emballeret i en indeslutning, der kan forsegles, eller den kryokonserverede hudækvivalent er tilvejebragt i en dyrkningsbeholder, der er emballeret i en pose.

**12.** Fremgangsmåde ifølge krav 1, endvidere omfattende at udforme den organotypisk dyrkede hudækvivalent før trin (a) ved:

at tilvejebringe en dyrkningsskål omfattende et cellekultursubstrat, der kan bevæges mellem definerede øvre og nedre positioner i dyrkningsskålen,

5 at udforme en hudækvivalent på cellekultursubstratet, hvor cellekultursubstratet er i den øvre position, og derefter sænke cellekultursubstratet til den nedre position til udførelse af trin (a) ifølge krav 1.

**13.** Fremgangsmåde ifølge krav 1, endvidere omfattende at udforme den organotypisk dyrkede hudækvivalent før trin (a) ved:

10 at tilvejebringe en dyrkningsskål omfattende en indsats, der kan bevæges mellem øvre og nedre positioner i dyrkningsskålen, hvor indsatsen har en plan bundflade, der er dannet af en porøs membran,

15 at udforme en dermal ækvivalent omfattende fibroblastceller på den porøse membran i indsatsen, hvor indsatsen er anbragt i den øvre position i dyrkningsskålen,

at dyrke fibroblastcellerne til dannelse af en dermal ækvivalent,

20 at påføre keratinocytter på den dermale ækvivalent, dyrke keratinocytterne i et dyrkningsmedium under betingelser, hvor keratinocytterne danner en hudækvivalent omfattende lagdelt epitel,

at fjerne dyrkningsmediet,

at sænke indsatsen til den nedre position til udførelse af trin (a) ifølge krav 1.

**14.** Emballeret kryokonserveret hudækvivalent, der er fremstillet ifølge krav 1, til anvendelse ved behandling af en patient, hvor den kryokonserverede hudækvivalent før påføring på patienten

a) optøs og ikke renses; eller

30 b) opvarmes, overføres aseptisk fra emballagen og bringes i kontakt med et absorberende medium omfattende en vævskompatibel opløsning for at muliggøre fjernelse af den kryobeskyttende opløsning ved diffusion.

**15.** Fremgangsmåde ifølge krav 1, hvor hudækvivalenten i trin (a) inkuberes i den kryobeskyttende opløsning i 10 til 60 minutter, fortrinsvis fra 20 til 30 minutter.



# DRAWINGS

Figure 1




Figure 2




Figure 3




Figure 4




Figure 5




Figure 6




Figure 7




Figure 8






Figure 9




Figure 10

