

(12) STANDARD PATENT
(19) AUSTRALIAN PATENT OFFICE

(11) Application No. **AU 2006237730 B2**

- (54) Title
Immobilisation of antigenic carbohydrates to support detection of pathogenic micro-organisms.
- (51) International Patent Classification(s)
G01N 33/569 (2006.01) **G01N 33/543** (2006.01)
- (21) Application No: **2006237730** (22) Date of Filing: **2006.04.24**
- (87) WIPO No: **WO06/112708**
- (30) Priority Data
- (31) Number (32) Date (33) Country
05075967.9 **2005.04.22** **EP**
- (43) Publication Date: **2006.10.26**
(44) Accepted Journal Date: **2012.03.15**
- (71) Applicant(s)
RnA Holding B.V.
- (72) Inventor(s)
Bokken, Gertruda Cornelia Antonia Maria;Gortemaker, Bertha Gerarda Maria;Bergwerff, Aldert Anthonie
- (74) Agent / Attorney
Freehills Patent & Trade Mark Attorneys, Level 43 101 Collins Street, Melbourne, VIC, 3000
- (56) Related Art
WO2000/036419 A1
Laferriere et al, (1997) Vaccine, vol. 15, pps 179-186
US5,116,962A

(19) World Intellectual Property Organization
International Bureau



(43) International Publication Date
26 October 2006 (26.10.2006)

PCT

(10) International Publication Number
WO 2006/112708 A1

(51) International Patent Classification:

G01N 33/569 (2006.01) **G01N 33/543** (2006.01)

Gerarda, Maria [NL/NL]; Karekiet 8, NL-3831 HM Leusden (NL).

(21) International Application Number:

PCT/NL2006/000218

(74) Agent: **VAN LOON, C.J.J.**; C/o Vereenigde, Johan De Wittlaan 7, NL-2517 JR The Hague (NL).

(22) International Filing Date: 24 April 2006 (24.04.2006)

(25) Filing Language: English

(26) Publication Language: English

(30) Priority Data:

05075967.9 22 April 2005 (22.04.2005) EP

(71) Applicant (for all designated States except US): **RNA HOLDING B.V.** [NL/NL]; Yalelaan 2, Kamer 307, NL-3584 CM Utrecht (NL).

(72) Inventors; and


(75) Inventors/Applicants (for US only): **BERGWERFF, Aldert, Anthonie** [NL/NL]; Dr. J.p. Heijelaan 1, NL-3818 GT Amersfoort (NL). **BOKKEN, Gertruda, Cornelia, Antonia, Maria** [NL/NL]; De Kersentuin 23, NL-6658 WJ Beneden-Leeuwen (NL). **GORTEMAKER, Bertha,**

(81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KM, KN, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, LY, MA, MD, MG, MK, MN, MW, MX, MZ, NA, NG, NI, NO, NZ, OM, PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SM, SY, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.

(84) Designated States (unless otherwise indicated, for every kind of regional protection available): ARIPO (BW, GH, GM, KE, LS, MW, MZ, NA, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IS, IT, LT, LU, LV, MC, NL, PL, PT,

[Continued on next page]

(54) Title: IMMOBILISATION OF ANTIGENIC CARBOHYDRATES TO SUPPORT DETECTION OF PATHOGENIC MICRO-ORGANISMS.



(57) Abstract: The invention relates to the field of chemistry and diagnosis, more in particular to diagnosis of current and/or past and/or symptomless infections or of a history of exposure to a gram-negative-bacterium (such as an enterobacteriaceae or a legionella). Even more in particular, the invention relates to the screening of animals or animal products for the presence of unwanted/undesired microorganisms. The invention further relates to a method for screening samples for the presence of antibodies directed against unwanted/undesired microorganisms and preferably such a method is performed with help of a biosensor. The invention also relates to a method for immobilising polysaccharides to solid surfaces. The invention furthermore provides solid surfaces with immobilised polysaccharides as well as applications of such surfaces.



RO, SE, SI, SK, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

Published:

— *with international search report*

Immobilisation of antigenic carbohydrates to support detection of pathogenic micro-organisms.

The invention relates to the field of chemistry and diagnosis, more in particular to diagnosis of current and/or past and/or symptomless infections or of a history of exposure to a gram-negative-bacterium (such as an enterobacteriaceae or a legionella). Even more in particular, the invention relates to the screening of animals or animal products for the presence of unwanted/undesired microorganisms. The invention further relates to a method for screening samples for the presence of antibodies directed against unwanted/undesired microorganisms and preferably such a method is performed with help of a biosensor. The invention also relates to a method for immobilising polysaccharides to solid surfaces. The invention furthermore provides solid surfaces with immobilised polysaccharides as well as applications of such surfaces.

The world is full of gram-negative bacteria, many of which are members of the family Enterobacteriaceae. Members of this family are found in the gastrointestinal tract of animals, but many are also free living in soil and water. Members of the family Enterobacteriaceae have very complex antigenic structures. Moreover, they comprise multiple antigens that are identified as K antigens, H antigens and O antigens. The K antigen is the acidic polysaccharide capsule. The capsule has many functions including evasion from the immune system of the infected host and adhesion to the epithelium of the host. The H antigen is located on the flagella.

The outer portion of the cell wall in gram-negative bacteria is chiefly composed of lipopolysaccharides (LPS). LPS is composed of lipid A which is buried in the outer membrane, a short carbohydrate core and optionally a chain of polysaccharides that is made up of repeating units. The O-antigens are located on the polysaccharide. Lipid A is the toxic constituent of the LPS. As cells lyse, LPS is released, leading to fever and complement consumption. It also interferes with coagulation and at high concentrations eventually leads to a state of shock.

As a non-limiting example one member of the enterobacteriaceae, salmonella, will be discussed in more detail. A large number of the subspecies of the genera of *Salmonella enterica* are important pathogenenic bacteria for humans and animals. Besides that animals go into a pathological episode, animals can be symptomless carriers of the bacteria. Contaminated animals can be a source of these pathogens threatening public health for example through the food that these animals produce. As many stakeholders consider the number of food-borne salmonella infections unacceptable, measures have to be taken to contain this pathogen in the food chain.

Salmonella is of major significance as a pathogenic microorganism in food-borne infections in humans, causing mild to severe clinical effects. In The Netherlands, 5% of all identified cases of gastroenteritis is salmonellosis (Edel et al., 1993; Hoogenboom Verdegaal et al., 1994). The average incidence of this infection is 450 cases per 100,000 person years at risk, which is similar to that in other industrialized countries (Berends et al., 1998). Despite the 2480 serotypes identified in the group of *S. enterica* up to 2001 (Popoff, 2001), only a small number have been involved in human infections (Grimont et al., 2000). *Salmonella typhimurium* plus *Salmonella enteritidis* represented >75% of all salmonella isolates from human sources sent to the Dutch National Salmonella Centre at the RIVM in 2002 (Van Pelt et al., 2003). This percentage consisted of 51% contributed by contact with chicken products (poultry 15%; eggs 36%)(Van Pelt et al., 2003).

Detection of immunoglobulins in the body fluids of organisms (serology) is a way to establish a history of exposure of animals and humans to infectious agents. A humoral response against salmonella antigens can be detected in chickens 1 week post-infection and persists for at least 10 weeks even if the bird is no longer culture-positive (Holt, 2000). The antigenic determinants of salmonella are, as described above, composed of somatic (O), flagellar (H) and surface (Vi) antigens (Holt, 2000). Variations in the composition of antigens correlate with different salmonella serotypes.

Typically, serology is faster than culture-typing of the disease-causative organism. Fast and specific detection of potential salmonella-positive

herds and flocks is of importance in order to take adequate measures in production processes. The detection of antibodies in serum and blood samples from food-producing animals reporting the presence of zoonotic pathogens is therefore of significance. Such information is then used as the input for risk-
5 assessment and rational slaughtering of potentially pathogen-contaminated animals in order to be able to increase food safety, but also to improve occupational hazards and to reduce spreading of the pathogens in the environment.

A number of serological tests have been developed for the detection of
10 invasive salmonella species. Among many such methods, agglutination and ELISA have most commonly been used (Barrow, 2000). Agglutination tests have been used successfully to eradicate *Salmonella pullorum* from poultry flocks. However, the approach is cumbersome, laborious and not suitable for large-scale screening programs according to modern standards. Several ELISA procedures,
15 which are considered relatively cheap and fast, have therefore been developed to detect anti-*S. enteritidis* and *S. typhimurium* antigen responses in poultry sera (Barrow et al., 1996; Thorns et al., 1996; de Vries et al., 1998; Barrow, 2000; Yamane et al., 2000).

The use of biosensors also promises to be useful, cheap and rapid in
20 this area of analysis. In addition, the technique is able to detect multiple analytes of any biomolecular type in a single run. A biosensor is defined as an analytical device consisting of (i) a re-usable immobilized biological ligand that 'senses' the analyte, and (ii) a physical transducer, which translates this phenomenon into an electronic signal.

25 The surface plasmon resonance (SPR) phenomenon was first recognized in the early 1960s (Kretschmann and Raether, 1968) and the first SPR biosensors were introduced in the 1980s (Liedberg et al., 1983). It took until the late 1980s and early 1990s before the first commercially available SPR-based biosensor equipment was released on the market. Initially, this type of biosensor
30 attracted the interest of pharmaceutical companies as a secondary tool for both selective and sensitive *in vitro* screening of promising novel pharmaceutical products from combinatorial libraries. It proved to be a valuable alternative for

2006237730 05 Jul 2011

classic approaches such as ELISA procedures. Moreover, it offers real-time measurement of the binding event in contrast to end-point determinations. The benefits of this analytical approach have also been recognized by many other life science disciplines, including food sciences (Ivnitski et al., 1999; Medina, 1997). So far, only a few publications on SPR biosensing have addressed the detection of pathogenic microorganisms, for example the use of immobilized *Escherichia coli* O157:H7 cells to screen the performance of anti- *E. coli* O157:H7 antibodies (Medina et al., 1997), and the use of these antibodies to detect *E. coli* O157:H7 cells (Fratamico et al., 1997).

In Jongerius-Gortemaker et al. (2002) a study to the suitability of an SPR optical biosensor to detect antibodies in serum and blood indicating a humoral reaction to invasion with *Salmonella* serotypes *enteritidis* and *typhimurium* was initiated. In this study, use was made of immobilised flagellar antigen fusion proteins. After thorough analysis it is concluded that the sensitivity and/or the robustness of this system was not sufficient and in particular not for high-throughput screening of for example poultry at the slaughter line in an abattoir, processing animals at the rate of several thousands per hour.

An aspect of the present invention is to provide for a method that has an improved sensitivity and/or an improved robustness. This goal has been reached by developing a carrier with immobilised somatic or so-called O-antigens. As described, the O-antigens are located on the lipopolysaccharides and the composition of the polysaccharide varies and corresponds with the serovar of the salmonella (sub)species. Every serotype can, amongst others, be described by a number of O-antigens and are typically coded with a number, such as O4, O6 or O12. The O-antigens can be found as repeating units on the polysaccharide part of the LPS. The length of the polysaccharide also varies and can be between zero (rough LPS) and more than 50 repeating units (smooth LPS).

Within the *Salmonella enterica* family, different serogroups can be distinguished; each group comprises at least one specific O-antigen. The salmonella serovars of importance in chicken and pigs are listed with their O-antigen profile in Table 1. In Denmark, Germany, Greece and The Netherlands, 39.5% of all salmonella-positive pigs sampled at the abattoir were determined as *S. typhimurium*. Dependent of country, other important isolates from pigs were *S. derby* (17.1%), *S. infantis* (8.0%), *S. panama* (5.1%), *S. ohio* (4.9%), *S. london* (4.4%), *S. livingstone* (3.1%), *S. virchow* (2.7%), *S. bredeny* (2.1%), *S. mbandaka* (1.1%), *S. Brandenburg* (1.0%), *S. goldcoast* (0.8%).

2006237730 05 Jul 2011

In case of chickens, 14% of the chickens were salmonella-positive at flock level in 2002 in The Netherlands. The predominant serovar was in that case *S. paratyphi B var. java*. At the retail level a comparable percentage (13.4%) was found in the Netherlands. The most frequent salmonella serovars isolated from broilers in 14 EU member states were *S. paratyphi B var. java* (24.7%), *S. enteritidis* (13.6%), *S. infantis* (8.0%), *S. virchow* (6.7%), *S. livingstone* (5.7%), *S. mbandaka* (5.5%), *S. typhimurium* (5.3%), *S. senftenberg* (5.0%), *S. hadar* (3.7%). *S. paratyphi B var. java* is dominating, but this is fully attributable to The Netherlands.

Reference to any prior art in the specification is not, and should not be taken as, an acknowledgment or any form of suggestion that this prior art forms part of the common general knowledge in Australia or any other jurisdiction or that this prior art could reasonably be expected to be ascertained, understood and regarded as relevant by a person skilled in the art.

As used herein, except where the context requires otherwise the term 'comprise' and variations of the term, such as 'comprising', 'comprises' and 'comprised', are not intended to exclude other additives, components, integers or steps.

Table 1. Some salmonella serovars considered as important zoonotic agents in broilers and in pigs listed with their O-antigen profiles (Popoff, 2001)

Salmonella serovar	Chicken (C)/ pigs (P)	O-antigen profile	serogroup
<i>Brandenburg</i>	P	4, [5], 12	B
<i>Bredeney</i>	P	<u>1</u> , 4, 12, <u>27</u>	B
<i>Derby</i>	P	<u>1</u> ^a , 4, [5] ^b , 12	B
<i>Enteritidis</i>	C	<u>1</u> , 9, 12	D ₁
<i>Goldcoast</i>	C/P	6, 8	C ₂
<i>Infantis</i>	C/P	6, 7, <u>14</u>	C ₁
<i>Livingstone</i>	P	6, 7, <u>14</u>	C ₁
<i>London</i>	P	3, 10, [15]	E ₁
<i>Mbandaka</i>	P	6, 7, 14	C ₁
<i>Meleagridis</i>	P	3, 10, [<u>15</u>], [<u>15</u> , <u>34</u>] ^c	E ₁
<i>Ohio</i>	P	6, 7, 14	C ₁
<i>Panama</i>	P	<u>1</u> , 9, 12	D ₁
<i>Paratyphi B var.</i>	C	<u>1</u> , 4, [5], 12	B
<i>Java</i>			
<i>Typhimurium</i>	C/P	<u>1</u> , 4, [5], 12	B
<i>Virchow</i>	P	6, 7, 14	C ₁

^a O antigen determined by phage conversion is indicated by underlining

^b O antigens which may be present or absent are indicated in square brackets

5 ^c lysogenized by phage ϵ_{15} [15] and by phage ϵ_{34} [15, 34]

In a first embodiment, the invention provides a method for immobilisation of a polysaccharide on a carrier, comprising contacting said polysaccharide with an oxidising agent and a polymer comprising at least two amine and/or amide groups to obtain a polysaccharide-polymer complex and coupling said

5 polysaccharide-polymer complex to said carrier. The polymer can be any polymer that contains at least two amine and/or amide groups. Said at least two amine and/or amide groups preferably cross-link said polymer to said polysaccharide and said carrier. To allow for more efficient coupling it is preferred that said polymer comprises at least 4 and more preferably at least 7 amine or amide
10 groups. The polymer comprises at least 10 building blocks. Building blocks of a polymer share characteristic reactive groups that enable elongation of the polymer. A preferred building block is an amino acid or a functional part, derivative and/or analogue thereof. In a preferred embodiment said polymer comprises a protein. A protein comprises at least one polypeptide chain
15 comprising at least 10 amino acids or functional equivalent thereof. A protein contains at least constituents having free amine and/or amide groups. In the context of the invention the protein can also be a multimer comprising at least two polypeptide chains that are covalently or non-covalently linked to each other. The protein may comprise modifications such as those common to biological
20 systems such as post-translational glycosylation. The protein may also be artificially modified or provided with a further group as long as it has the mentioned amine and/or amide groups available.

In a preferred embodiment said polysaccharide is derived from a gram-negative bacterium. The sensitivity of such a prepared carrier is much improved
25 when the lipopolysaccharide (O-antigen) before the immobilisation on the carrier is oxidised in the presence of a polymer comprising at least two amine and/or amide groups, preferably a protein. Although we do not wish to be bound by any theory, it is currently thought that the aldehyde groups that result from the oxidation of the polysaccharides are capable of reacting with the amino groups of
30 the protein to form a substituted imine (Schiff-base binding). Upon injection over (an activated) carrier (for example a sensorchip) the available aldehyde groups react with hydrazide to form hydrazon. The following reduction stabilises not

only the covalent binding between the carrier (for example a carrier comprising dextran) and the polysaccharide but also the imine binding between protein and polysaccharide. As will be explained in more detail in the experimental part, polysaccharides (O antigens) of different salmonella sera types have been
5 immobilised on a carrier. The prepared carriers were subsequently subjected to an SPR-analysis with standard sera. The obtained serological response was used as an indicator for success of the method. When coupling reactions were performed without the oxidation step no or almost no significant response of reference sera could be detected.

10 Preferably, the immobilisation/coupling of the polysaccharide-protein complex to a carrier is such that high sensitivity and/or robustness is obtained. Whereas flagellar antigens denature and lose their antigenicity towards serum antibodies while the sensor chip has to be regenerated for a next analysis cycle with relatively harsh solvents, the somatic antigens are found rather stable
15 towards these regeneration solvents. In fact, the loss of immobilized O-antigen activity is believed to be primarily associated with degradation of the solid surface, namely gradual loss of dextran layer attached to the goldfilm, to which the antigens are bound. The method according to the invention results in a carrier that is more robust compared to a carrier of the prior art.

20 Preferably, the invention provides a method for immobilisation of a polysaccharide on a carrier, comprising contacting said polysaccharide with an oxidising agent and a protein to obtain a polysaccharide-protein complex and coupling said polysaccharide-protein complex to said carrier, wherein said
25 polysaccharide is derived from a gram-negative bacterium and even more preferably wherein said polysaccharide is derived from an enterobacteriaceae. Yet even more preferably, said polysaccharide is derived from a gram-negative bacterium that is a human or veterinary or plant pathogen. Examples of such polysaccharides are polysaccharides derived from a salmonella (sub)species.
30 Other examples are polysaccharides derived from *Eschericia coli* species (for example E.coli O157) and the bacterial species outlined in Table 2.

Table 2. Examples of LPS-containing bacteria pathogenic to human and/or animals.

Bacterial species	Mainly found in	affecting
<i>Campylobacter coli</i>	Swine	Humans
<i>Campylobacter jejuni</i>	Avian species, dogs	Humans
<i>Campylobacter lari</i>	Seagull	Humans
<i>Escherichia coli</i> O157	Ruminants	Humans
<i>Legionella pneumophila</i>	Water	Humans
<i>Listeria monocytogenes</i>	As food process- contaminant	Humans
<i>Salmonella choleraesuis</i>	Swine	Swine
<i>Salmonella enteritidis</i>	Avian species, swine	Humans
<i>Salmonella gallinarum</i>	Avian species	Chickens
<i>Salmonella goldcoast</i>	Swine	Humans
<i>Salmonella infantis</i>	Chickens	Humans
<i>Salmonella livingstone</i>	Swine	Humans
<i>Salmonella meleagridis</i>	Swine	Humans
<i>Salmonella pollorum</i>	Avian species	Chickens
<i>Salmonella typhimurium</i>	Avian	Humans, horses
<i>Streptococcus suis</i>	Swine	Swine, humans
<i>Vibrio cholerae</i> (non-O1)	Aquatic animals	Humans
<i>Vibrio parahaemolyticus</i>	Aquatic animals	Humans
<i>Vibrio vulnificus</i>	Aquatic animals	Humans
<i>Yersinia enterocolitica</i>	Swine	Humans

As will be explained in more detail later, a carrier comprising an
 5 immobilised polysaccharide (O-antigen) is particularly useful in the diagnosis of
 the mentioned LPS-containing bacteria.

The term "polysaccharide" is intended to mean an entity comprising
 two or more glycoside linked monosaccharide units and embraces, amongst
 others, an oligosaccharide (2-10 residues) as well as a polysaccharide (more than

10 monosaccharides). The linking may result in linear or branched polysaccharide. In a preferred embodiment, the invention provides a method for immobilisation of a polysaccharide on a carrier, comprising contacting said polysaccharide with an oxidising agent and a protein to obtain a polysaccharide-protein complex and coupling said polysaccharide-protein complex to said carrier, wherein said polysaccharide is a lipopolysaccharide (LPS), i.e. a polysaccharide comprising lipid A. It is clear to a skilled person that the used (lipo)polysaccharide must comprise at least one antigenic structure and one group available/suitable for providing a linkage between the protein and the polysaccharide. More details in respect of the last item will be provided later on. Hence, as long as the (lipo)polysaccharide comprises an antigenic structure and a group suitable for providing a linkage between the protein and the polysaccharide an immobilization method of the invention may be used to obtain a sensitive and/or robust carrier.

The LPS is expressed at the cellular exterior and is part of the bacterial cellular wall. The expression of LPS is not under direct genetic control, so that LPS is a pool of different molecules with varying composition of the lipid A part in terms of the attached aliphatic chain. Moreover, the bacterial cell may synthesize rough LPS with no or a short carbohydrate chain, or smooth LPS with a mature carbohydrate chain existing of more than 50 repeating units expressing its antigenicity. In addition to this heterogeneity, within a single molecule LPS, several O-antigen entities, which are distinctively numbered, may be expressed. An O-antigen profile is, however, per definition unique for a salmonella serogroup. A complete serotyping of a salmonella also includes the H-antigens as well as the Vi-antigens.

LPS may be obtained by a variety of methods and the experimental part describes in more detail the use of a trichloric acid extraction (optionally followed by ethanol extraction and dialysis) according to Staub (1965) for this purpose. Other examples of suitable extraction methods are described by Wilkons (1996) and include, but are not restricted to, extractions with diethylene glycol, dimethyl sulphoxide, NaCl-diethyl ether (1:2 (v/v)), NaCl-butan-1-ol (1:1 (v/v)),

aqueous EDTA, NaCl-sodium citrate, aqueous phenol or aqueous phenol-chloroform petroleum.

The purity of the obtained/used LPS batch is considered not to be extremely critical. It is experienced that the LPS does not have to be completely free of contaminants. The specific coupling reaction provides a certain degree of selectivity. Moreover, as described in the experimental part, the used/obtained LPS (preferably an LPS batch) is optimised in respect of the amount of protein necessary for an optimal response. It is clear to a skilled person that the LPS preferably comprises not much rough LPS. The preferred LPS batch essentially comprises smooth LPS.

Although we do not wish to be bound by any theory it is currently thought that the presence of a 2-keto-3-deoxy-octonic acid (KDO) and/or a glycerol-mannoheptose (Hep) and/or a GlcNAc in the core of the LPS molecule is needed for a covalent coupling.

Although a lot of different bacteria are employed by the term gram-negative bacteria it is believed that LPS from all these bacteria are suitable for use in the presently claimed invention as long as the LPS comprises at least one constituent with non-conjugated or de-conjugated vicinal hydroxy groups, preferably in the core region of the LPS molecule. In salmonella, most likely candidate constituents are KDO and Hep and GlcNAc residues. The presence or absence of such a KDO and/or Hep and/or GlcNAc group is indirectly genetically determined. Although the genetic information necessary to construct the species-, serotype or even strain specific monosaccharides is present in the corresponding organism, it depends on the growth circumstances whether said LPS contains aldehyde-convertible monosaccharides in the core region.

There are of course also other sources of LPS available, such as buying it commercially.

In a preferred embodiment, the invention provides a method for immobilisation of a polysaccharide on a carrier, comprising contacting said polysaccharide with an oxidising agent and a protein to obtain a polysaccharide-protein complex and coupling said polysaccharide-protein complex to said carrier, wherein said protein is a protein (for example a serum protein) with a certain

amount of (primary) amines. Preferably, at least some of these amines are not sterical hindered and/or are not participating in non-covalent bindings, such as H-H bridges or dipole-dipole interactions and/or are not protonated. Such a protein preferably does not have or hardly have, any immunogenic properties and hence cross-reacting antibodies directed to the used protein are avoided as much as possible. Examples of suitable proteins are haemoglobin (Hb), ovalbumin (Ob), myoglobin (Mb) and serum albumin (SA). The biosensor response of different standard sera on immobilised LPS oxidized in the presence of Hb or Ob or Mb or SA were determined. Serum albumin, myoglobin and haemoglobin gave the most promising results. In a preferred embodiment the protein is haemoglobin or myoglobin.

The necessary protein is obtained commercially or by (overexpressing) in a suitable expression system or by isolating it from a suitable source. Haemoglobin has for example been obtained by isolating it from blood. Preferably the used protein batches are as pure as possible, thereby circumventing as much cross-reactions as possible. It is however experienced that small amounts of contamination are allowed without jeopardising the sensitivity and/or robustness of the obtained carriers.

The ratio (lipo)polysaccharide versus protein depends, amongst others, on the used protein. Experiments with Hb have shown that concentrations between 15 and 50 % (m/m) have resulted in satisfactory results. When bovine serum albumin is used much lower ratios, between 0.7 and 7% (m/m), are used. Some examples: the optimal Hb concentration for *S. livingstone* LPS is around 50% (m/m) and for *S. enteritidis* LPS the optimal Hb concentration is 15% (m/m).

The isolated LPS preparations are preferably oxidised in the presence of a protein facilitated by an oxidising agent. In a preferred embodiment the invention provides a method for immobilisation of a polysaccharide on a carrier, comprising contacting said polysaccharide with an oxidising agent and a protein to obtain a polysaccharide-protein complex and coupling said polysaccharide-protein complex to said carrier, wherein said oxidising agent is capable of oxidising vicinal diols. Even more preferably, the oxidising agent preferably

oxidises vicinal diols at least under controlled condition. Oxidation of vicinal diols is preferred as this warrants reliable coupling of vicinal diol containing polysaccharide to the matrix. In a preferred embodiment of the invention the polysaccharides to be coupled to the matrix contain an antigen that is to be

5 recognised by a member of a binding pair. To be recognizable it is preferred that the antigen is left unchanged at least in the majority of the polysaccharides that are being coupled to the carrier. This requires a balance between the level of oxidation required to obtain efficient coupling to the matrix and availability of the antigen for association with the member of the binding pair. The latter

10 requires that the antigen is left essentially unaffected by the oxidation at least in an amount sufficient to be usable in a diagnostic setting. Oxidation of vicinal diols according to the present invention warrants the availability of sufficient antigen in recognizable form while at the same time allowing efficient coupling of the polysaccharide to the carrier. In a preferred embodiment said oxidising agent

15 comprises (sodium) *m*-periodate. Other periodates such as potassium periodate or other salts thereof are also suitable periodates of the present invention. Periodate oxidation is very suited for enabling preferential oxidation of vicinal diols according to the present invention. Oxidation of predominantly vicinal diols in a polysaccharide of the invention can typically be achieved by incubating said

20 polysaccharide with said periodate at a concentration of between 1 and 10 mM periodate. Other parameters of the reaction influence both the speed and the type of reaction predominantly performed. One example is incubation time. When applying very short incubation times higher than 10 mM periodate can be used. Periodate preferably oxidises vicinal diols, particularly of the more susceptible

25 vicinal diols in the side chains of the polysaccharide. Thus as long as so-called 'mild' reaction conditions are chosen, preferably vicinal diols will be oxidised. When conditions are chosen that also allow other oxidation reactions to occur more often (for instance because of depletion of the vicinal diol substrate), the antigen present in the polysaccharide will be affected significantly. Thus for the

30 present invention a periodate oxidation is said to be mild when the mentioned preferred concentrations are used and when at least 20% and preferably at least 50%, more preferably at least 70% and most preferably about 90% of the antigen

is intact after oxidation. Availability or intactness of the antigen is preferably measured by means of an ELISA assay using a standardized antibody. Again we do not wish to be bound by any theory but it is currently thought that periodate will induce an oxidative disruption of linkages between vicinal diols on especially carbohydrate moieties, as in *e.g.* mannose, to yield aldehyde functionalities. This reaction is typically performed in buffers at a pH range between 4.5 and 5.5 in the dark using a (preferably) freshly prepared 1-100 mM sodium meta-periodate in 0.1 M sodium acetate. Preferably the reaction is performed at a concentration of between 1 and 10 mM metaperiodate. The oxidation is performed in the presence of a protein in the ranges as discussed above. The bis-aldehyde compounds, like the oxidised monosaccharide constituents in the polysaccharide chain of LPS, may react with any amino group in a protein and may form a Schiff-base linkage resulting in a substituted imine. When one or both of the vicinal hydroxyl groups is condensed in a covalent sugar linkage, the hydroxyl function is lost and no oxidation occurs. This is the case in many branched and/or linearly linked oligo- and polysaccharides. In the case of salmonella LPS, the inner core structure carries in most cases an oxidisable Gal, GlcNAc, Hep and/or KDO, but non-reducing Hep and KDO constituents are most susceptible for oxidation, in particular at very mild oxidation conditions at concentrations less than 6 mM meta periodate. Because the core region is a rather conserved part of LPS from different Enterobacteriaceae, (lipo)polysaccharides of members of the Enterobacteriaceae may be applied in a method of the invention.

Periodate will also oxidise, when present, certain aminoethanol derivatives such as the hydroxylysine residues in collagen, as well as methionine (to its sulfoxide) and certain thiols (usually to disulfides). In addition, N-terminal serine and threonine residues of peptides and proteins can be selectively oxidized by periodate to aldehyde groups. These reactions, however, usually occur at a slower rate than oxidation of vicinal diols and the presence of such group does not substantially interfere with a method according to the invention.

The invention also provides a method for immobilisation of a polysaccharide on a carrier, comprising contacting said polysaccharide with an oxidising agent and a protein to obtain a polysaccharide-protein complex and

coupling said polysaccharide-protein complex to said carrier, further comprising a step which results in ending/stopping the oxidation process, for example by desalting of said polysaccharide-protein complex. This is for example accomplished with help of a NAP-5 column. However the person skilled in the art is aware that many other methods exist which have the same effect, for example adding a reductor or an easily oxidisable molecule such as glycerol. Preferably, the way of stopping the oxidation is such that at the same time a buffer change is accomplished, for example HPLC, FPLC, dialysis, ion-exchangers, gel electrophoresis or ultrafiltration.

For storage purposes, the production of evaporated aliquots, after addition of protein, is also described within the experimental part. This results in the presence of a large stock of reproducible material.

The invention therefore further comprises the obtained intermediate, i.e. the preparation of in the presence of protein oxidised polysaccharide, optionally desalted and optionally evaporated.

Preferably, the used carrier is made of an inert, non-hydrophobic material and the binding of the LPS-protein complex to said carrier is covalent. Even more preferred such a carrier has a low protein binding or low biomolecular binding. Examples are a carrier of glass or silica or of a non-hydrophobic plastic. In a preferred embodiment said carrier is in the form of a microsphere or bead. Several types of microsphere or beads are available to the person skilled in the art. In a preferred embodiment said microsphere or bead comprises polystyrene. Microsphere or beads are particularly preferred because they can be provided with different antigens using a method of the invention. Microsphere or beads with different antigens can be accordingly coded with a different colour. Testing a sample for the presence of an antibody against an antigen can be done using a collection of the mentioned microsphere or beads. Binding of the antibody to a particular type of antigen can now be detected easily by the colour code of the microsphere or bead bound. Binding of the antibody can be detected in various ways. For instance, microsphere or beads containing bound antibody can be extracted from the sample and measured using a further antibody specific for the

constant region of the antibody. On the other hand, sample can be directly analysed, i.e. in the absence of further manipulations by labelling the bound antibody and simultaneously detecting colour of the antibody and the colour of the microsphere or bead. Various methods for simultaneous detection of two or
5 more colours are available to the person skilled in the art. In the present invention, a colour is defined as any type of electromagnetic radiation that can be detected, be it a typical colour revealed, for instance, by reflection of light, to light emitted as a result of fluorescence or phosphorescence.

The invention thus further provides a collection of at least two
10 microsphere or beads wherein at least two of said at least two microsphere or beads each comprise a different antigen of the present invention. In a preferred embodiment said antigen comprises O-antigen of Salmonella. In a particularly preferred embodiment said antigen is linked to said microsphere or beads carrier using a method of the invention. Thus preferably at least one of said microsphere
15 or beads comprises a polysaccharide coating linked to a polysaccharide comprising an antigen to be detected linked to each other via a polymer comprising at least two amine and/or amide groups, preferably a protein of the invention, wherein said linkage polymer (protein) is linked to said polysaccharide comprising said antigen, via an amine and/or amide group on said polymer and a
20 periodate oxidised vicinal diol on said polysaccharide comprising said antigen.

In a preferred embodiment the invention provides a method for immobilisation of a polysaccharide on a carrier, comprising contacting said polysaccharide with an oxidising agent and a protein to obtain a polysaccharide-
25 protein complex and coupling said polysaccharide-protein complex to said carrier, further comprising activating the surface of said carrier. In an even more preferred embodiment, said carrier comprises a glass surface coated with gold and even more preferred said carrier is modified with a carboxyl donor. A surface can be activated. Carboxylic acid (COOH) groups (further referred to as carboxyl
30 groups) are needed on this surface. Preferably these COOH groups are provided by a stable homogeneous layer of molecules, which may have been modified for this purpose. These surfaces may exist of, but are not limited to, carboxylic acid-

modified polysaccharides, alkanes or alkenes, such as polyethylene, attached to e.g. gold, polystyrene or silicon surfaces. Preferably the carrier comprises a polysaccharide that acts as a carboxyl donor. More preferably a carboxymethylated dextran layer wherein said polysaccharide modified carrier, preferably comprising a dextran layer is activated with 1-ethyl-3-(3-dimethylaminopropyl)carbodiimide hydrochloride and N-hydroxysuccinimide. The activation is preferably followed by preparation with carbodiazine. In the next step the polysaccharide-protein complex is added to the activated dextran layer. The reactive aldehyde functionalities react spontaneously with the hydrazide to hydrazones, which are then reduced to stabilise the covalent bonds.

Prior to routinely use, the performance of chip-conjugated LPS to bind anti-Enterobacterium (for example salmonella) antibodies is assessed using reference polyclonal agglutination sera.

Depending on the analytical/diagnostic question asked it is decided whether one or for example at least two different serogroup-representing carbohydrates, preferably 4 different serogroup-representing carbohydrates are used. In case one is interested in knowing which particular serogroup is present, multiple (the amount of which is different on the particular question asked and on the used apparatus) different serogroup-representing carbohydrates are used and if one just wants to know whether for example an animal is or has been infected by a particular serogroup, a single serogroup-representing carbohydrate may be oxidised in the presence of a protein and immobilised on a carrier. The use of at least two different serogroup-representing carbohydrates results in a carrier that can be used in a multi-serogroup analysis. More preferably at least three and even more preferred at least more than three (for example four) different polysaccharides are used. These polysaccharides may be oxidised in the presence of one type of protein or in the presence of different types of protein. The skilled person is capable of making any sensible combination. For example, to be able to detect more than 90% of all salmonella infections serogroups B, C and D in chicken and serogroups B, C, D and E in pigs should be represented.

Using one type of serogroup-representing carbohydrates is extremely useful if one is interested in the question whether or not a certain type of bacterium is or was present. Using multiple different serogroup-representing carbohydrates is for example useful if one wants to determine whether an animal is or was infected by any gram-negative bacteria (for example enterobacteriaceae).

In a preferred embodiment the invention provides a method for immobilisation of a polysaccharide on a carrier, comprising contacting said polysaccharide with an oxidising agent and a protein to obtain a polysaccharide-protein complex and coupling said polysaccharide-protein complex to said carrier, wherein said carrier is a biosensor chip. Such a biosensor chip is commercially available (for example that produced by Biacore) and hence no further information will be provided.

In another embodiment the invention provides a carrier obtained by the method according as described above or a carrier comprising an immobilised polysaccharide-protein complex on its surface. In one embodiment of the invention a carrier of the invention comprises a polysaccharide coating that is linked to a further polysaccharide coating via reductive amination, wherein said further polysaccharide coating comprises a protein coupled to said further polysaccharide coating through oxidation of vicinal diols on said further polysaccharide - protein. In a preferred embodiment said reductive amination is achieved. In a preferred embodiment the invention provides a carrier comprising a polysaccharide coating that is coupled to polysaccharide

In yet another embodiment the invention provides biosensor comprising a carrier according to the invention. Whether the carrier is obtained by a method according to the invention can for example be determined by extracting the polysaccharides from said carrier and determining whether covalently linked protein is present. As already discussed above the carrier may

also comprise different immobilised polysaccharides (for example O-antigens) possibly in combinations with different types of protein. However, also one type of protein may be used in the oxidation of different polysaccharides.

Whether a carrier and/or chip of the invention is employed can for example be determined with help of MALDI-MS possibly in combination with proteolytic digestion. Such an analysis provides information with respect to the used protein and polysaccharide. With help of acidic hydrolysis the polysaccharide-protein complexes are released from the carrier. Such an obtained mixture is then subjected to LC-MS/MS analysis before and after proteolytic hydrolysis. The obtained complex may also be subjected to a monosaccharide analysis, for example GC-MS following methanolysis and/or Smith degradation, from which it is determined which type of LPS is used. This information is furthermore used to determine whether KDO, Hep or other sugars have been oxidised.

15

A carrier of the invention may be used in different detection systems, for example optical, thermal, acoustic, amperometric, magnetic or chemical and a carrier of the invention may be used in any biomolecular interaction assay (BIA) or any affinity assay (AA). As a non-limiting example, the use of optical detection via Surface Plasmon Resonance is described in more detail.

20

The invention provides a Surface Plasmon Resonance detection system comprising a biosensor as described above. The gold layer in the sensor chip creates the physical conditions required for Surface Plasmon Resonance (SPR). The principle of SPR will be described in the context of Biacore instruments. They incorporate the SPR phenomenon to monitor biomolecular interactions in 'real-time'. At an interface between two transparent media of different refractive index such as glass and water, light coming from the side of higher refractive index is partly reflected and partly refracted. Above a certain critical angle of incidence no light is refracted across the interface and total internal reflection (TIR) occurs at the metal film-liquid interface. This is where light travels through an optically dense medium such as glass, and is reflected back through that medium at the interface with a less optically dense medium such as buffer.

30

Although the incident light is totally reflected, the electromagnetic field component, termed the evanescent wave, penetrates a distance on the order of one wavelength into the less optically dense medium. The evanescent wave is generated at the interface between a glass prism (high refractive index) and a layer of buffer (lower refractive index). If the interface between the media of higher and lower refractive indices is coated with a thin metal film (a fraction of the light wavelength), then the propagation of the evanescent wave will interact with the electrons on the metal layer. Metals contain electron clouds at their surface, which can couple with incident light at certain angles. These electrons are also known as plasmons, and the passage of the evanescent wave through the metal layer causes the plasmons to resonate, forming a quantum mechanical wave known as a surface plasmon. Therefore, when surface plasmon resonance occurs, energy from the incident light is lost to the metal film resulting in a decrease in the reflected light intensity. The resonance phenomenon only occurs at an acutely defined angle of the incident light. This angle is dependent on the refractive index of the medium close to the metal-film surface. Changes in the refractive index of the buffer solution (e.g. an increase in surface concentration of solutes), to a distance of about 300 nm from the metal film surface will therefore alter the resonance angle. Continuous monitoring of this resonance angle allows the quantitation of changes in refractive index of the buffer solution close to the metal-film surface. In 'real-time' Biacore, the metal film properties, wavelength, and refractive index of the glass (denser medium) are all kept constant, and as a result SPR can be used to monitor the refractive index of the aqueous layer immediately adjacent to the metal (gold) layer. In the Biacore system the chip is composed of glass, has 4 channels and the associated gold layer is covered with a layer of dextran chemically modified to facilitate immobilisation of ligands such as antibodies or antigens. Any changes in mass that occur due to binding of the analyte with the immobilised antibody on the sensor chip will cause a change in SPR angle, which is monitored in 'real-time' and quantified as a sensorgram. A mass change of approximately 1 kRU (1,000 RU) corresponds to a mass change in surface protein concentration of 1 ng/mm². Typical responses for surface binding of proteins are of the order of 0.1-20 kRU.

There is no need to label molecules with fluorescent or radioactive tags- so avoiding the possibility that labels may compromise activity and moreover no difficult or expensive chemistry is necessary for labelling.

The obtained carriers can be used in different types of analysis, such as
5 bacteriology (direct assay) or serology (indirect assay).

An example of a serological assay is a method for determining the presence of an antibody directed to an antigen of a gram-negative bacteria in a sample, comprising contacting said sample with a carrier or a biosensor as described above and determining whether the carrier has bound any antibody
10 (Figure 1).

Such a method is for example very suitable for determining the presence of an antibody directed against an O antigen and thus it is indirectly established whether an infection is present or whether a recent infection has occurred. Such a method is for example used to screen slaughter animals for salmonella or to
15 screen animals for salmonella before they are exported abroad. Moreover, the method is also applied to samples obtained from living (for example, farm or zoo) animals.

Examples of samples that can be used in such a method are tissue sample, body fluid, secretes or excretes and more detailed examples are blood,
20 blood derived samples, tissue, meat juice, milk, egg, fluids from an eye, saliva or faeces. As already outlined the samples can be obtained from dead as well as living animals.

A method according to the invention is not limited to a certain immunoglobulin (sub)type but can in principle be every (iso)type immunoglobulin
25 such as (s)IgA₁, (s)IgA₂, IgD, IgG₁, IgG₂, IgG₃, IgG₄, IgM, IgY. Moreover, it may also be any other antigen-binding material. Preferably, such an antigen-binding material is a biomarker of a (history) of infection.

Such a serological assay is for example directed to one particular serogroup-representing carbohydrate or to different (i.e. multi analyte)
30 serogroup-representing carbohydrates and hence such a method is for example used to determine the presence or absence of a certain salmonella (sub)type.

An example of a bacteriological assay is a method for determining the presence of a gram-negative bacterium in a sample, comprising contacting said sample with a predetermined amount of antibodies directed against an antigen of said bacterium and determining the amount of antibodies not bound to said bacterium with a carrier or a biosensor as described above.

Preferably the antigen is a serogroup-representing carbohydrate.

This method optionally further comprises the removal of non-bound antibodies from contacted sample and predetermined amount of antibodies by for example washing or immuno-magnetic separation procedures, centrifugation or filtering.

For this type of analysis every type of sample can be used, such as animal feed, manure, feathers, soil, water for consumption or sewage water, meat, orange juice, chocolate, skin, vegetables etc. Animal samples may be obtained from living as well as dead animals.

In this bacteriological assay a single type of antibody as well as a mixture of at least two different types of antibodies (directed against different antigens, for example two different serogroup-representing carbohydrates) is used.

Preferably, such serological and bacteriological assays are performed such that the binding to said carrier or said biosensor is determined by Plasmon Surface Resonance or fluorescent microsphere or bead counter.

The source of the samples is as already outlined above unlimited and may for example be obtained from a human or an animal. Examples of suitable animals are (race) horses, pigs, poultry (for example chicken, turkey, quail, duck, and goose), ruminants (for example calf or cow, goat, sheep). The animals may be farm animals, zoo animals as well as free living animals. Moreover, samples from these animals may be obtained from living as well as dead animals.

In yet another embodiment the invention provides a method for determining the presence of a gram-negative bacterium in a sample comprising

- contacting said sample with target bacteria-specific, bacteriophages and allowing the bacteriophages to infect said sample
- removing non-bound and/or non-invading bacteriophages resulting in a bacteriophage infected sample
- 5 - bringing the bacteriophage infected sample into contact with an indicator organism susceptible for the used bacteriophages
- incubate during at least one bacteriophage multiplication cycle
- recover the bacteriophages to obtain a bacteriophage-containing sample
- analyse said bacteriophage-containing sample with a carrier or a biosensor
- 10 according as described above.

Most analytical methods require prior enrichment and growth in specific media to detect bacteria, including salmonella. Usually sample preparation is very time-consuming relatively to the total analysis time. It generally takes 3 to 5 days before the presence of e.g. salmonella can be
15 confirmed. In many situations, this time for analysis is unacceptable and hinders trade and indirectly threatens community health.

The objective of this part of the invention is development of a fast (preferably within 24 h) and/or cost-effective and/or specific and/or sensitive diagnostic method for the determination of the presence of micro organisms. For
20 this reason, the development of a biomolecular interaction assay (BIA) which exploits the ability of genus- and/or serovar-specific bacteriophages to multiply in their 'victim' bacteria, is aimed. An increment in number of the target pathogen-specific phage(s) indicates not only the presence of the target organism but is also a (semi-)quantitative measure for the content of target bacteria in the tested
25 sample.

A schematic overview of the proposed BIA method is depicted in Figure 2. A particulate sample is homogenised for example using a Stomacher. Liquid samples are mixed by vigorous shaking. Analyte cells are then extracted or enriched by any suitable method and may comprise (a combination of) selective
30 growth, centrifugation, filtration and/or immuno-magnetic separation (IMS). Enriched cells are fortified with target bacteria-specific bacteriophages and incubated for a few minutes while mixing. Before the multiplication cycle of the

bacteriophage is complete, cells are washed to remove as complete as possible any non-bound and non-invading bacteriophages. Following the multiplication cycle of the bacteriophage, the sample is brought in contact with an indicator organism susceptible, i.e. in a life phase that is sensitive for bacteriophage penetration and intracellular multiplication, for the used bacteriophage, preferably at the highest possible concentration (for example concentrated overnight culture). The bacteriophage-bacterium suspension is incubated for at least one bacteriophage multiplication cycle. The phage-infected suspension is then centrifuged or filtered to precipitate/ remove cellular material and to recover multiplied bacteriophages. The bacteriophage-containing sample is injected over an LPS-conjugated biosensor chip (according to the invention) to retain these particles in the detector for the generation of analyte-specific biosensor response.

To gain as much time as possible the indicator organism can be kept as a continuous culture in the lab and has a cell density of usually 10^9 CFU/ml.

Such a suspension may be concentrated to 10^{10} CFU/ml, as higher cell densities will increase sensitivity of the proposed method.

This method can be used to determine a single type of serovar but to detect multiple serovars in one run, a mix of different bacteriophages and a mixture of possibly different indicator bacteria may have to be applied.

Target bacteria-specific bacteriophages are described in the prior art and examples are provided in the experimental part, for example anti-*Salmonella enteritidis* bacteriophages.

Phages have been described to attach to LPS, including the phage described in the experimental part for salmonella detection. Suitable carriers/chips are carriers/chips with LPS or with immobilised bacterial surface molecules (thus including membrane proteins and other biomolecules or a combination thereof). Use of LPS of cell membrane material will circumvent the generation of poly- or monoclonal antibodies. If attachment of the phages to bacterial biomolecules (LPS) is not satisfactory in the BIA, biosensor chip-immobilised anti-phage antibodies may have to be used in a successful BIA to capture bacteriophages from the probed sample.

The invention furthermore provides a kit with components suitable for use in any of the described applications. Depending on the customer's demand, such a kit comprises a ready-for use carrier/chip obtained by a method according to the invention. When the customer wants to prepare the carrier himself, the kit will at least comprise (lipo)polysaccharide fortified/enriched with protein (for example haemoglobin or serum albumin) in a predetermined amount, an amount of oxidizing agent (for example periodate), suitable buffers. Optionally, such a kit comprises means for desalting, for example a desalting column. When the customer wants to mix (lipo)polysaccharide and protein himself these components are delivered separately together with an instructions manual. Optionally, such a kit may furthermore comprise positive and/or negative reference sera, a sample dilution buffer and any necessary instruction manual.

The methods as described above are particularly suitable for screening samples on a large-scale basis. In one of the earlier (slow) settings 96 samples were checked within 33 minutes. In a large-scale setting with relative slow biosensor equipment 15.000 samples were screened within 3 months. This number could have been much higher but unfortunately one of the slaughterhouses stopped participating.

The invention will be explained in more detail in the following description, which is not limiting the invention.

Examples

Example 1

5 **Materials and methods**

1.1 Materials

1.1.1 Chemicals

10 Amine coupling kits, consisting of *N*-hydroxysuccinimide (NHS), 1-ethyl-3-(3-dimethylaminopropyl)carbodiimide hydrochloride (EDC) and ethanolamine hydrochloride – sodium hydroxide pH 8.5 and the running buffer (HBS-EP), containing 10 mM HEPES, 150 mM sodium hydrochloride, 3 mM EDTA and 0.005% (v/v) surfactant P20 at pH 7.4, were bought from Biacore AB (Uppsala, Sweden), which also supplied ready-to-use 10 mM glycine and 50 mM sodium hydroxide. Ethanol, ethylene glycol, sodium chloride, sodium hydroxide and trichloroacetic acid (TCA) were purchased from Merck (Darmstadt, Germany). Carboxymethylated-dextran sodium salt, sodium cyanoborohydride and carbohydrazide were obtained from Fluka Chemie GmbH (Buchs, Switzerland).
20 CHAPS (Plus one) was delivered by Pharmacia Biotech (Uppsala, Sweden). Sodium acetate trihydrate and acetic acid were supplied by J.T. Baker (Deventer, The Netherlands). Guanidine hydrochloride was obtained from Calbiochem (San Diego, CA, U.S.A.). *Porcine* haemoglobin (Hb) and myoglobin (Mb), chicken ovalbumin (Ob; 98% grade V), *bovine* serum albumin (BSA; 96% Fraction V),
25 sodium periodate, Tween-20, Tween-80 and Triton X-100 were acquired from Sigma Chemical Company (St. Louis, MO, U.S.A.). Water was obtained from of a Milli Q water purification system (Millipore, Bedford, MA, U.S.A.).

1.1.2 Materials

NAP-5 columns (0.5 ml; Sephadex G-25) were purchased from Amersham
30 Biosciences (Roosendaal, The Netherlands) and were used as described by the producer. CM5 biosensor chips were bought from Biacore AB. Dialysis bag

(Spectra/Por) with a cut-off of 1 kDa was obtained from Spectrum Laboratories Inc. (Rancho Dominguez, CA, U.S.A.)

1.1.3 Anti-salmonella antisera

The following salmonella monovalent 'O' somatic *lapine* antisera were used: anti-O4, anti-O5, anti-O6,7, anti-O8, anti-O9, anti-O10, anti-O12, O Poly E (anti-O3, anti-O10, anti-O15, anti-O19, anti-O34). In addition, salmonella polyvalent 'O' somatic (Poly A-S) *lapine* antisera (anti-O2, anti-O3, anti-O4, anti-O5, anti-O6,7, anti-O8, anti-O9, anti-O10, anti-O11, anti-O12, anti-O13, anti-O15, anti-O16, anti-O17, anti-O18, anti-O19, anti-O20, anti-O21, anti-O22, anti-O23, anti-O28, anti-O30, anti-O34, anti-O35, anti-O38, anti-O40, anti-O41) was used as well. The sera were purchased from Pro-Lab diagnostics (Salmonella Reference Section of the Central Veterinary Laboratory, Weybridge, U.K.). Serogroup specific *murine* anti-B (anti-O4, O5 en O27), anti-C (anti-O7, O8), anti-D (anti O9, Vi) and anti-E (anti-O3, O19) monoclonal antibodies were bought from SIFIN (Berlin, Germany').

Sera were diluted 1:20 (v/v) in HBS-EP containing 1.0 M sodium chloride, 1% (m/v) carboxymethylated dextran and 0.05% (v/v) Tween 80, except anti-O5 serum was diluted 1:200 (v/v) and the anti-serogroup specific preparations were diluted 1:100 (v/v) in the same solvent.

1.1.4 Reference avian and porcine sera

All reference sera were obtained from the Dutch Animal Health Service (Deventer, The Netherlands). The obtained avian reference sera were reactive with *Salmonella enteritidis* (serogroup D₁), *S. typhimurium* (serogroup B), *S. pullorum/gallinarum* (serogroup D₁) and *S. infantis* (serogroup C₁), and were further referred to as C-Se, C-St, C-Spg and C-Si, respectively. These chicken sera were originally prepared for ELISA analyses as positive references. In addition, specific pathogen-free chicken serum (further referred to as C-SPF) was purchased as a negative control reference sample. These sera were reconstituted from freeze-dried material by addition of water at a volume indicated by the manufacturer. C-Se, C-Spg and C-Si were diluted 1:200 (v/v) in HBS-EP containing 1.0 M sodium chloride, 1.0% (m/v) carboxymethylated dextran and

0.05% (v/v) Tween-80, whereas C-SPF and C-St were diluted 1:50 (v/v) in the same solution. Likewise, *porcine* sera from animals challenged with *S. typhimurium* and *S. livingstone* (serogroup C₁) were referenced as P-St and P-Sl, respectively. In addition, *Actinobacillus pleuropneumoniae* serotype 2-
5 reacting *porcine* serum used as control in a complement fixation test, was exploited as negative control for porcine serum in the salmonella biosensor assay. The *porcine* sera were diluted 1:20 (v/v) in HBS-EP containing 1.0 M sodium chloride, 1% (m/v) carboxymethylated dextran and 0.05% (v/v) Tween 80 as end concentrations.

10 1.1.5 Salmonella stock

The bacteria *Salmonella goldcoast* (Sg; serogroup C₂), *S. livingstone* (Sl) and *S. melaegridis* (Sm; serogroup E₁) were obtained from an in-house collection, while *S. enteritidis* #23 phage type Pt4 (Se), and *S. typhimurium* X-193 phage type 507 (St) were kind gifts of F. G. van Zijderveld (Animal Sciences Group,
15 Lelystad, The Netherlands). The bacteria were grown in overnight cultures in Nutrient Broth #2 (Oxoid, Basingstroke, U.K.). Stocks of salmonella strains were morphologically and biochemically confirmed as salmonella and also verified for the presence of the correct, expected O-antigens by an agglutination reaction of the cells with specific standard anti O-antigen anti-sera (Pro-Lab diagnostics) as
20 indicated in Table 3 on a glass plate. After addition of a half of the original volume with glycerol (Merck), stocks were stored in portions at -80°C.

Table 3. *O*-antigen verification of the salmonella serovars used for LPS production. The expected reaction of anti sera used for verification by agglutination, is given

Salmonella Serovar	Agglutination sera								
	α -O4 ^a	α -O5	α -O6,7	α -O8	α -O9	α -O10	α -O12	poly A-S	Poly E
<i>S. enteritidis</i> (O9,O12) ^b	-	-	-	-	+	-	+	+	-
<i>S. goldcoast</i> (O6,8)	-	-	+	+	-	-	-	+	-
<i>S. livingstone</i> (O6,7)	-	-	+	-	-	-	-	+	-
<i>S. meleagridis</i> (O3,O10)	-	-	-	-	-	+	-	+	+
<i>S. typhimurium</i> (O4,O5,O12)	+	+	-	-	-	-	+	+	-

^a α -O4: antibodies reacting with antigen structure coded with O4; in a similar way the antibodies against O5, O6,7, O9, O10 and O12 are indicated.

^b O antigens specific for salmonella serovar is indicated in brackets

1.2 Methods

1.2.1 Extraction of LPS

Overnight cultures of salmonella were prepared by applying 100 µl from their
5 corresponding stocks on each of the 120 plates containing brain heart infusion
agar (BHIa, Oxoid). The presence of the expected salmonella serovar was
confirmed through conventional selective growth, bio- and immunochemical
classification, whenever new stock suspensions were produced. The bacteria
were harvested from the surface of the plates into 1 ml 9 g/l NaCl (saline)
10 solution per agar plate using a trigalski spatula. Each plate was washed twice
with 2 ml saline solution. Bacteria were collected in six centrifugation tubes.
Each tube was complemented with 100 ml saline and mixed before
centrifugation at 10,000 g and 4°C for 15 min and supernatant was discarded.
This centrifugation step was repeated twice by suspending cells in 75 ml saline
15 wash solution per tube each run. While kept on ice, pelleted bacteria were
suspended in water at a volume ratio, which was a 5-fold to the weight of the
bacteria. An equivolume of 0.250 M (Se) or 0.500 M (Sg, Sl, Sm and St) TCA
was added to give end concentrations of 0.12 M and 0.25 M, respectively,
followed by continuous stirring at 4°C for 3 h. A lipopolysaccharide (LPS)-
20 containing supernatant was then acquired at 20,000 g and 4°C for 30 min. The
pH of the supernatant was adjusted to pH 6.5 with 5 M sodium hydroxide and
when nearing the aimed pH with 0.10 M sodium hydroxide. The final volume
of the LPS-containing solution was determined prior to storage at -18°C for
30 min. The solution was diluted with a double volume of freezing cold
25 absolute ethanol from a -18°C storage place, and incubation was continued
overnight at -4°C without stirring in a closed, in house-build device with
circulating cold ethylene glycol/water (1:4, v/v). An LPS-containing pellet was
obtained after centrifugation at 20,000 g and -4°C for 30 min. The particulate
material was suspended in a volume of 0.5 ml water per gram original
30 bacterial mass weighed at the start of extraction process. The suspension was

dialyzed in a 1-kDa dialysis bag against water at 4°C for two days with regular intermittent refreshment of the water. The bag content was centrifuged at 20,000 *g* and at 4°C for 30 min, and the supernatant was lyophilized. The lyophilisate was weighed to establish the recovery of LPS. LPS was

5 reconstituted in water to make up an end concentration of 5 mg/ml.

Dependent of type of LPS and batch (see also section 1.2.2), a volume of 1 mg/ml *porcine* haemoglobin (Hb) was added to a concentration as indicated in the text. Each batch was portioned into 0.5-mg LPS fractions, which were dried using a vacuum evaporator and then stored at 5-8°C.

10

1.2.2 Optimal haemoglobin content

Protein was added to an LPS preparation prior to its chemical modification and immobilization to a sensor chip to acquire high coating levels and high serum responsive antigens. The optimum Hb content in each LPS batch was

15 established by comparison of the responses of immobilized LPS that was fortified with Hb at different levels, using a panel of positive and negative reference sera.

1.2.3 Oxidation of LPS

20 A portion of 0.5 mg haemoglobin-fortified LPS was dissolved in 500 µl 100 mM sodium acetate pH 5.5. Following the addition of 20 µl 50 mM sodium periodate, the solution was incubated for 40 min on ice protected from light. The oxidation of LPS was quenched and the solution was desalted by passing 500 µl of the reaction mixture through an NAP-5 cartridge with a gravity-

25 controlled flow. Modified LPS was eluted with 1 ml 10 mM sodium acetate, pH 4.0. Prior to use, the cartridge was conditioned thrice with 3 ml 10 mM sodium acetate, pH 4.0.

1.2.4 Immobilization of LPS

To immobilize the antigens to a sensor chip, the following handlings were conducted at a flow rate of 5 μ l/min in a Biacore 3000 instrument controlled by Biacore 3000 Control Software (version 3.1.1; Biacore). Immobilization of
5 oxidized LPS was achieved by execution of the aldehyde-coupling procedure described in BIAapplications Handbook, version AB (1998). Briefly, the dextran layer at the biosensor chip CM5 was activated with a 7-min pulse of a mixture of EDC/NHS available from the amine-coupling kit. The activation was immediately followed by injection of 5 mM aqueous carbonyldiimidazole for
10 7 min as well.

Deactivation of the excess of reactive groups was then accomplished with a pulse of 1 M ethanolamine for 7 min. Prior to immobilisation of the antigen, LPS was diluted in sodium acetate pH 4.0 in a ratio dependent of the salmonella serovar (see text) and immobilised for 32 min. The linkage
15 between dextran-matrix and antigen was then stabilized by injection of 100 mM sodium cyanoborohydride solved in 10 mM sodium acetate at pH 4 at a flow rate of 2 μ l/min for 20 min. A relative response indicative for a successful LPS immobilisation procedure is 2 kRU for a 62.5 μ g/ml LPS solution containing 15% (m/m) protein, and 9 kRU for a 250 μ g/ml LPS
20 solution containing 50% (m/m) protein.

1.2.5 SPR biosensor assay

Optical SPR biosensor assays were performed on a Biacore 3000 SPR biosensor platform controlled by the same software as described above. Prior to
25 injection, sera were diluted in HBS-EP buffer containing 1.0% (m/v) carboxymethylated-dextran sodium salt, 1.0 M sodium chloride and 0.05% (m/v) Tween 80 at a ratio of 1:50 (v/v) or otherwise as indicated in the text. The mixtures were incubated for at least 2 min at ambient temperature. Pig sera were injected for 2 min at 40 μ l/min, whereas bird sera were injected for
30 2 min at 5 μ l/min or 20 μ l/min as indicated.

Regeneration of the chip to recover the antigenic activity of the sensor surface was achieved with a 15-s pulse of 6 mM glycine at pH 2, containing 6 M guanidine hydrochloride, 0.1% (m/v) CHAPS, and 0.1% (v/v) of each Tween-20, Tween-80 and Triton X-100. This was followed with a second regeneration
5 step with the running HBS-EP buffer enriched with 0.05% (m/v) CHAPS (end concentration) for 12 s at 100 μ l/min.

1.2.6 Monosaccharide analysis

Trimethylsilylated (methyl ester) methyl glycosides were prepared from the
10 glycan samples by methanolysis (1.0 M methanolic HCl, 24 h, 85°C) followed by re-*N*-acetylation and trimethylsilylation, and then analyzed by gas chromatography/mass spectrometry as described [Kamerling JP, Vliegthart JFG (1989)]. The quantitative analysis was carried out by gas chromatography on a capillary EC-1 column (30 m $\times$ 0.32 mm, Alltech) using a
15 Chrompack CP 9002 gas chromatograph operated with a temperature program from 140°C to 240°C at 4°C/min, and flame-ionization detection. The identification of the monosaccharide derivatives was confirmed by gas chromatography/mass spectrometry on a Fisons Instruments GC 8060/MD 800 system (Interscience) equipped with an AT-1 column (30 m $\times$ 0.25 mm,
20 Alltech).

Results

LPS isolation

For the production of LPS, yields of bacterial cells and of LPS were compared
5 for agar plate culture and growth of salmonella in broth (Table 4). For
laboratory technical reasons, it was decided to harvest bacteria from agar
plates, rather than isolation of the cells from culture flasks. The results of the
isolation of LPS from Se, Sg, Sl, Sm and St are summarized in Tables 5 to 9,
respectively. The standardized isolation of well-defined LPS is determinative
10 for a successful and robust serological assay. To secure assay performance,
batch-to-batch differences should be kept to a minimum. For this reason,
several batches of LPS extracted from each Se, Sg, Sl, Sm and St were
produced. The recovery of LPS largely depended on the final TCA
concentration in the mixture during extraction of LPS (*cf.* Table 6 and Table
15 8), although this relationship was not completely clear for the extraction of
LPS from St (Table 9). Indeed, no accurate optimal TCA concentration could
be determined for each LPS type through the testing of a broad range of TCA
concentrations. Here, optimal TCA would yield highest LPS amounts, and
give highest specific serological and lowest aspecific biosensor responses. In
20 this study, the TCA concentration chosen as 'optimal' for LPS extraction from
the different salmonella serotypes was based on the final LPS yields after
dialysis, and were 0.12 M, 0.25 M, 0.25 M, 0.25 M and 0.25 M as end
concentrations for Se, Sg, Sl, Sm and St, respectively.

Table 4. Recovery of LPS from *S. enteritidis* cells grown either as a suspension in a bioreactor containing so-called nutrient broth#2 (broth) or on BHI agar plates (agar). LPS was isolated using indicated TCA end concentrations. The yield of LPS relative to the amount of isolated cells is indicated in the last column.

LPS Batch code	Culture method	TCA (M)	bacteria yield (g)	recovered LPS (mg)	LPS yield (%, m/m)
Se01	Broth	0.25	3,9	16	0,42
Se02	Broth	0.25	5,0	21	0,41
Se03*	Agar	0.25	14	0,2	0,00
Se04a	Broth	0.25	3,4	0,4	0,01
Se04b	broth	0.5	3,4	2	0,06
Se05**	agar	0.5	7,6	2,4	0,03
Se06a	agar	0.5	8,8	3,9	0,04
Se06b	agar	0.25	7,6	9,5	0,12
Se06c	agar	0.125	8,9	13	0,14
Se07a	agar	0.1	9,1	12	0,13
Se07b	agar	0.05	9,2	2,9	0,03
Se07c	agar	0.025	9,3	4	0,04
Se2003.1	agar	0.125	29	48	0,16
Se2003.2	agar	0.125	17	25	0,15
Se2003.4	agar	0.125	13	5,1	0,04
Se2005.1	agar	0.25	15	18	0,13

* pH of TCA-containing mixture is outlying

** some material was lost during sample work up process.

Table 5. Recovery of *Salmonella enteritidis* cells grown on BHla plates. LPS was isolated using 0.12 M TCA end concentration (cf. Table 4). Rec, recovered.

LPS Batch (Se)	Total bacterial yield (g)	Number of BHla plates	bacteria per plate (g)	Rec. LPS/cells (% m/m)
Se2003.1	29.09	98	0.29	0.16
Se2003.2	17.27	60	0.28	0.15
Se2003.4	13.00	40	0.32	0.04
Se2005.1	44.5	120	0.37	0.12

Table 6. Recovery of *Salmonella goldcoast* cells grown on BHla plates. LPS was extracted using a TCA end concentration as indicated. Rec, recovered.

LPS Batch code	Total bacterial yield (g)	Number of BHla plates	bacteria per plate (g)	TCA ^a (M)	Rec. LPS/cells (% m/m)
Sg2003.1	40.39	120	0.34	0.075	0.01
Sg2003.2	35.09	120	0.29	0.25	0.41
Sg2003.3	12.89	40	0.32	0.25	0.36
Sg2005.1	46.99	120	0.39	0.25	0.30 ^b

^a end concentration TCA in extraction mixture.

^b approximately a third of the production was lost during work-up.

Table 7. Recovery of *Salmonella livingstone* cells grown on BHIa plates. LPS was extracted using 0.250 M TCA end concentration. Rec, recovered.

LPS Batch code	Total bacterial yield (g)	Number of BHIa plates	bacteria per plate (g)	Rec. LPS/cells (% m/m)
Sl2003.1	32.89	120	0.27	0.52
Sl2003.2	13.63	40	0.34	0.51
Sl2005.1	47.40	120	0.40	0.64

5 **Table 8.** Recovery of *Salmonella meleagridis* cells grown on BHIa plates. LPS was isolated using a TCA end concentration as indicated. Rec, recovered.

LPS Batch code	Total bacterial yield (g)	Number of BHIa plates	Bacteria per plate (g)	TCA ^a (M)	Rec. LPS/cells (% m/m)
Sm2003.1a ^b	9.10	30	0.30	0.250	0.32
Sm2003.1b ^b	10.13	30	0.34	0.125	0.19
Sm2003.1c ^b	10.00	30	0.33	0.075	0.06
Sm2003.2 ^b	40.42	120	0.33	0.075	0.02
Sm2003.3	37.28	138	0.27	0.250	0.47

^a end concentration TCA in extraction mixture.

^b batches Sm2003.1 to 2003.2 were combined to a single batch called Sm2003.1

Table 9. Recovery of *Salmonella typhimurium* cells grown on BHIa plates. LPS was isolated using 0.250 M TCA end concentration. Rec, recovered.

Batch code	Total bacterial yield (g)	Number of BHIa plates ^a	bacteria per plate (g)	TCA ^a (M)	Rec. LPS/cells (% m/m)
St2003.1	28.51	69	0.41	0.125	0.18
St2003.2 ^b	39.21	120	0.33	0.250	0.06
St2003.3 ^b	16.6	56	0.30	0.125	0.09
St2003.4a ^b	18.15	60	0.30	0.250	0.18
St2003.4b ^b	19.2	60	0.32	0.125	0.06
St2005.1	46.80	120	0.39	0.250	0.19

^a end concentration TCA in extraction mixture;

^b batches St2003.2 to St2003.4b were combined to a single batch called

5 St2003.2

The monosaccharide composition of isolated LPS preparations were analyzed to reveal the consistency of the isolation and purification procedure for LPS from different salmonella growths. It must be noted that analyses were

10 performed on LPS preparations that were ready for oxidation and for that reason fortified with Hb at levels that were determined most optimal for the LPS batch tested (see below). For this purpose, GC-FID and GC-MS analyses were carried out after methanolysis of the Hb-fortified LPS preparations (Table 10 through Table 14). These results show that Hb does not contribute

15 to a significant amount of carbohydrates in the final LPS preparation.

Analysis of BHIa, showed the presence of exclusively galactose (Gal) and glucose (Glc). The content of these monosaccharides was 5.6 µg/mg dried BHIa. Analyses of the salmonella LPS preparations, demonstrated the occurrence of Gal, Glc, *N*-acetyl glucosamine (GlcNAc), glycerol-manno-heptose

20 (Hep), 2-keto-3-deoxy-octonic acid (KDO), mannose (Man) and rhamnose (Rha);

6-deoxy-mannose) in accordance with their carbohydrate structures. Their relative occurrence was expressed as a molar ratio relative to 1.0 Man (as part of the PS region) or relative to 3.0 Hep (as part of the core region). It should, however, be noted that the core region contains 2 or 3 Hep residues.

5 Furthermore, GlcNAc can originate from either GlcNAc as in the repeating unit of SI LPS, or from glucosamine (GlcN), which occurs as disaccharide in the lipid A moiety as backbone for the attached lipids. Gal occurs in the core region, which is conserved in all *Salmonella enterica* serovars, and in the PS region of Se, Sg, St and Sm as well. In these cases, the molar ratio of Gal is
10 expected to be in excess of 1.0 Man, except for Sg in which each repeating unit contains 2 Man residues. The monosaccharide analyses did not include the detection of O-acetylated, phosphoryl-ethanolaminated or phosphorylated constituents, nor that of abequose (Abe; 3,6-dideoxy-xylohexose) or tyvelose (Tyv; 3,6-dideoxy-arabinose), which occur in the polysaccharide and core
15 regions of the isolated LPS types as well.

Analysis of Se LPS, showed the occurrence of Gal, Man and Rha at a molar ratio of 1.4, 1.0 and 1.2, respectively, in batch Se2003.1, whereas this ratio was 1.1, 1.0 and 0.9, respectively, in batch Se2003.4 (Table 5). This ratio is in good compliance with the composition of a repeating unit as [Tyv-]Man-Rha-Gal,
20 except the Rha ratio was significantly too high in batch Se2003.1. The carbohydrate content calculated on the basis of determined monosaccharides, was significant higher in batch Se2003.4, namely 241 µg compared 123 µg of batch Se2003.1. Considering the occurrence of 2 GlcN residues in the lipid A and a single GlcNAc residue in the core region and a single Man residue in
25 each repeating unit, the number of repeating units was estimated 19 and 20 in batches Se2003.1 and Se2003.4, respectively.

Monosaccharide analysis of oxidized Se2003.1 clearly demonstrates significant differences with the non-oxidized identical batch (Table 10). In contrast to the two GlcN residues, it is expected that the non-reducing, terminal GlcNAc

residue is for the greater part oxidized. Alditol derivatives were not detected by the monosaccharide analysis applied and a corresponding amount of GlcNAc-ol was not determined. As Gal and Man in the repeating unit are not susceptible towards periodate oxidation, the molar ratio of Man in the oxidized
5 batch is normalized to that of Man in the non-oxidized batch. It should be noted that both Gal residues in the core region are susceptible towards oxidation and thus the total Gal ratio is affected. Inspection of the molecular structure of Se LPS, suggests that in addition to terminal GlcNAc and core Gal, terminal KDOII-KDOIII disaccharide, and conjugated HepI and terminal
10 HepIII are susceptible to periodate oxidation as well. Indeed, the molar ratios of these monosaccharide residues suggest the loss of one Hep residue and approximately 1.6 KDO residues. It should be noted that KDOII may not be completely oxidized when this residue is conjugated with a phosphoryl-ethanolamine group.

Table 10. Monosaccharide analysis of *S. enteritidis* LPS and of oxidized *S. enteritidis* LPS. LPS was fortified with Hb at 15% (m/m). Molar ratios were determined on the basis of two GlcN and one GlcNAc residues (detected as three GlcNAc residues) present in the core and lipid A regions (referred to as CORE) and on the basis one Man residues in the repeating unit (referred to as UNIT). Normalized GlcNAc and Man residues are indicated by underlining. Carbohydrate content was determined in 0.5 mg LPS preparations, except monosaccharide analysis was performed on 125 mg oxidized material.

Monosaccharide	Molar ratio					
	Batch Se2003.1		Batch Se2003.1		Batch Se2003.4	
	CORE	UNIT	CORE	UNIT	CORE	UNIT
Gal	29.7	1.4	28.7	1.3	26.3	1.1
Glc	6.3	0.3	7.5	0.4	6.9	0.3
GlcNAc	<u>3.0</u>	+	2.2	+	<u>3.0</u>	+
Hep	2.8	+	1.9	+	2.0	+
KDO	3.2	+	1.4	+	2.3	+
Man	21.6	<u>1.0</u>	<u>21.6</u>	<u>1.0</u>	23.5	<u>1.0</u>
Rha	24.8	1.2	24.4	1.1	20.9	0.9
Carbohydrate content (µg) ^a	123.0		- . ^b		241	
Nr of repeating units	19		- -		20	

^a does not include the contribution of Tyv residues;

^b Amount of LPS-containing material analysed was not accurately determined.

Compared to Se LPS, monosaccharide analysis of Sg LPS suggests that LPS structure were smaller as the number of repeating units was significant lower (Table 11). The relative contribution of core Gal to PS Gal is for that reason larger and total molar ratio is found 1.5. In a similar way, the molar ratio for

Glc is found at 1.3 (batch Sg2003.2) and 1.5 (batch Sg2003.3), whereas the molar ratio for Rha fits with the expected structure. Batch Sg2003.2 seems however to contain less terminal HepIII and terminal KDOIII and could therefore offer less possibility for immobilization to the sensor chip.

- 5 **Table 11.** Monosaccharide analysis of LPS isolated from *S. goldcoast* fortified with Hb at 50% (m/m). Molar ratios were determined on the basis of two GlcN and one GlcNAc residues (detected as three GlcNAc residues) present in the core and lipid A regions (referred to as CORE) and on the basis two Man residues in the repeating unit (referred to as UNIT). Normalized GlcNAc and
- 10 Man residues are indicated by underlining. Carbohydrate content was determined in 0.5 mg LPS preparations.

Monosaccharide	Molar ratio			
	Batch Sg2003.2		Batch Sg2003.3	
	CORE	UNIT	CORE	UNIT
Gal	15.3	1.5	14.8	1.5
Glc	13.5	1.3	14.3	1.5
GlcNAc	<u>3.0</u>	+	<u>3.0</u>	+
Hep	2.9	+	2.3	+
KDO	3.0	+	2.8	+
Man	20.5	<u>2.0</u>	19.2	<u>2.0</u>
Rha	10.6	1.0	10.2	1.0
Carbohydrate content (μg) ^a	170		199	
Nr of repeating units	9		8	

^a does not include the contribution of Abe residues.

- Normalisation of the number of core residues from the monosaccharide analysis results of Sl LPS (Table 12) was hampered by the occurrence of
- 15 GlcNAc in the repeating units of the PS. When the number of Hep residues

was set at 3.0, an unacceptable overestimation of the number of KDO residues arose. For that reason, the number of Hep was set at 2.0, but may need to be modified, so that the number of KDO is closer to 3.0. As the number of Man residues in each repeating unit is four, molar ratios were corrected for 4.0 Man
5 residues. On the basis of a molar ratio of 4:1 of Man/GlcNAc in the PS region, the number of repeating units was calculated on the basis of the remaining core GlcNAc and Lipid A GlcN residues. This calculation revealed that the number of repeating units in SI was also relatively small, namely 8 and 10 units in batch SI2003.1 and batch2003.2, respectively.

Table 12. Monosaccharide analysis of LPS isolated from *S. livingstone* fortified with Hb at 50% (m/m). Molar ratios were determined on the basis of two Hep residues present in the core (referred to as CORE) and on the basis of four Man residues in the repeating unit (referred to as UNIT). Normalized GlcNAc and Man residues are indicated by underlining. Carbohydrate content was determined in 0.5 mg LPS preparations. n.d., not detected.

Monosaccharide	Molar ratio			
	Batch SI2003.1		Batch SI2003.2	
	CORE	UNIT	CORE	UNIT
Gal	3.5	0.4	5.1	0.4
Glc	17.4	1.7	23.5	1.6
GlcNAc	14.0	1.4	18.7	1.3
Hep	<u>2.0</u>	0.2	<u>2.0</u>	0.14
KDO	2.7	0.3	2.7	0.2
Man	40	<u>4.0</u>	57	<u>4.0</u>
Rha	n.d.	n.d.	n.d.	n.d.
Carbohydrate content (µg)	212		239	
Nr of repeating units	8		10	

Table 13. Monosaccharide analysis of *S. meleagridis* LPS. Batches Sm2003.1 and Sm2003.3 were fortified with 50% (m/m) Hb. Molar ratios were determined on the basis of two GlcN and one GlcNAc residues (detected as three GlcNAc residues) present in the core and lipid A regions (referred to as CORE) and on the basis one Man residues in the repeating unit (referred to as UNIT). Normalized GlcNAc and Man residues are indicated by underlining. Carbohydrate content was determined in 0.5 mg LPS preparations.

Monosaccharide	Molar ratio			
	Batch Sm2003.1		Batch Sm2003.3	
	CORE	UNIT	CORE	UNIT
Gal	20.4	1.5	22.4	1.4
Glc	7.9	0.6	4.8	0.3
GlcNAc	<u>3.0</u>	+	<u>3.0</u>	+
Hep	2.8	+	2.9	+
KDO	2.8	+	2.9	+
Man	14.2	<u>1.0</u>	15.4	<u>1.0</u>
Rha	16.1	1.1	17.6	1.2
Carbohydrate content (μ g) ^a	170		220	
Nr of repeating units	12		13	

^a does not include the contribution of *O*-acetyl groups, which may be attached to the repeating Gal residues.

- 10 Monosaccharide analysis of Sm LPS (Table 13) showed a completely different composition as that for Sl LPS in accordance with its molecular structure containing Man-Rha-Gal repeating units. As described above, the molar ratio for Gal is more than the expected 1.0 in the repeating unit partly by the contribution of Gal residues in the core region.

Likewise, equimolar ratios are expected for Glc, Man, Rha and Gal as these residues form a repeating unit in St LPS (Table 14). It should be noted here that abequose is also part of the repeating unit, but is not in the analysis applied. Batch St2003.2, however, contains less oxidizable Hep and KDO,
5 which may affect the efficacy of the immobilization of LPS from this preparation. At the other hand, Batch St2003.2 contains much more carbohydrate than batch St2003.1, namely 249 µg relative to 154 µg in 0.5 mg LPS, respectively.

Table 14. Monosaccharide analysis of *S. typhimurium* LPS. Batches St2003.1 and St2003.2 were fortified with Hb at 15% (m/m) and 25% (m/m), respectively. Molar ratios were determined on the basis of two GlcN and one GlcNAc residues (detected as three GlcNAc residues) present in the core and lipid A regions (referred to as CORE) and on the basis one Man residues in the repeating unit (referred to as UNIT). Normalized GlcNAc and Man residues are indicated by underlining. Carbohydrate content was determined in 0.5 mg LPS preparations.

Monosaccharide	Molar ratio			
	Batch St2003.1		Batch St2003.2	
	CORE	UNIT	CORE	UNIT
Gal	24.7	1.4	24.2	1.4
Glc	17.0	1.0	16.2	1.0
GlcNAc	<u>3.0</u>	+	<u>3.0</u>	+
Hep	2.8	+	2.6	+
KDO	2.7	+	2.4	+
Man	20.5	<u>1.0</u>	19.2	<u>1.0</u>
Rha	19.9	1.1	19.5	1.2
Carbohydrate content (µg) ^a	154		249	
Nr of repeating units	15		14	

^a does not include the contribution of (*O*-acetylated) Abe residues.

Protein-supported immobilization of LPS

Unexpectedly, intact LPS poorly coupled through its KDO-carboxylic acid function to EDC/NHS-activated carboxymethylated dextran, and no significant responses of reference sera were observed. To improve immobilization of LPS
5 to the sensor chip, LPS was oxidized using sodium periodate to create reactive aldehyde groups in its carbohydrate constituents, which would allow the so-called aldehyde coupling procedure, *i.e.* condensation of aldehyde with a hydrazide function into a hydrazone linkage followed by reduction to a hydrazide product. Without oxidation, LPS had indeed little potential to
10 immobilize to a surface of a CM5 chip coated with carbohydrazide (results not shown).

Coupling of oxidized LPS, however, gave disappointing reactivity with reference sera, probably as a result of insufficient immobilization of the
15 antigens. Commercially acquired phenol-extracted LPS, either intact or detoxified (*i.e.* cleavage of lipid A), from *S. enteritidis* gave low responses when immobilized after oxidation, namely 187 RU and 167 RU, respectively. It must be noted, however, that besides poor coupling, oxidation may have destroyed a part of the antigenic structures, which may give poor serological responses.
20 The degree of oxidation was investigated by monosaccharide analysis of Se LPS (Table 10). This analysis revealed that relative amounts of KDO, Hep and GlcNAc, which are constituent of the core region and not of the repeating antigenic units in the PS part, were significantly reduced compared to non-oxidized Se LPS. Importantly, monosaccharide residues part of the PS, and
25 thus antigenic structures, apparently remained intact under the mild oxidation conditions, which were applied.

Table 15. Biosensor responses following immobilization and responses of reference antisera flowed over chip surfaces, which were prepared with oxidized and non-oxidized Se or St LPS in the presence of 15% (m/m) *porcine* Hb.

Type of LPS	Periodate treatment	Biosensor response (RU)					
		Level of immobilization	α -O4 ^a	α -O5	α -O9	α -O12	O poly A-S
Se	No	6515	1	9	6	9	4
St	No	4402	2	156	1	5	0
Se	yes	4265	-21	-3	77	202	159
St	yes	5950	302	5005	-12	225	137

5 ^a anti-serum against indicated O antigen was tested

Table 16. Biosensor responses following immobilization and responses of reference antisera flowed over chip surfaces, which were prepared with oxidized St LPS (batch St2003.1) in the presence of 7% (m/m) of the indicated protein.

Protein added	Biosensor response (RU)											
	Level of immob. ^a	O4	O5	O9	O12	O poly E	O poly A-S	C-SPF	C-St	C-Se	C-Si	C-Spg
BSA ^b	344	n.d. ^c	n.d.	n.d.	n.d.	n.d.	n.d.	n.d.	n.d.	n.d.	n.d.	n.d.
BSA	7450	19	216	1	18	4	14	5	18	21	12	42
Hb	3410	319	3874	4	192	3	203	11	151	140	51	762
Mb	815	59	773	3	45	2	47	8	46	34	20	152
Ob	5540	76	926	10	56	12	56	12	47	45	25	177

5 ^a level of immobilization;

^b BSA was added to oxidized and desalted LPS;

^c immobilization of LPS was considered too low for further reference sera analysis.

- 10 To optimize the binding of LPS and consequently improve detection of binding antibodies from sera, oxidation of LPS was then executed in the presence of a protein to allow the formation of protein-LPS complexes through Schiff-base reactions between proteinaceous amines and aldehyde functions of LPS. Indeed, commercially available TCA extracted Se LPS, containing considerable
- 15 amounts of bacterial proteins, gave improved immobilization at 562 RU compared to phenol-extracted and ion-exchange chromatography-purified Se LPS at 208 RU.
- Oxidation of LPS was necessary, as mixtures containing non-oxidized LPS and protein show relatively high immobilization levels but insignificant specific
- 20 responses (Table 15). In addition, protein addition was only beneficial prior to oxidation of LPS, as addition of BSA to oxidized and desalted St LPS gave acceptable immobilization levels but no expected serological responses (Table

16). In a similar way, Hb yielded relatively high immobilization levels but no serological responses, as expected (results not shown).

For method improvement purposes, proteins with relatively high degree of
5 homology of their primary and secondary structure between homeothermic
vertebrate species and occurring in serology-suitable matrices were selected for
further investigations. For that reason, the performance of chicken ovalbumin,
porcine haemoglobin, *bovine* serum albumin or *porcine* myoglobin fortified (7%,
m/m) St LPS was compared (Table 16). Haemoglobin gave clearly best
10 improvement of immobilization levels together with best expected antigen-
antibody reactivity profile. In the presence of Hb, in particular, O12, poly O
A-S and C-St reference sera gave better responses.

This experiment was repeated with the addition of BSA, Hb and Mb at levels
15 indicated in Tables 17 to 20 using batches Se2005.1, Sg2005.1, Sl2005.1 and
St2005.1. This time, O4 and O5 bound to immobilized St LPS as expected
(Table 20). Evaluation of these results summarized in Tables 16 to 20 revealed
that when considering all expected responses simultaneously per LPS type, the
addition of Hb gave highest specific responses compared to the addition of BSA
20 and Mb. Furthermore, in most cases standard deviations that occurred with
Hb as supportive protein, were more favorable than those for the addition of
BSA and Mb. Haemoglobin was, therefore, selected for further
experimentation.

It was observed that the O poly A-S anti sera probably contains a low anti-
25 serogroup C₁ and C₂ titers, as in the case of testing immobilized Sg LPS
(Table 18) and Sl LPS (Table 19) relatively low responses are found.

Table 17. Biosensor responses in response units (RU) following immobilization and responses of reference antisera flowed over chip surfaces, which were prepared with Se LPS oxidized in the presence of 15% (m/m) of the indicated protein. Values were corrected for the C-SPF responses, which are listed as well. Standard deviations are indicated in brackets ($N = 5$, except for chicken and swine sera $N = 4$).

protein	Immobilization level ^a	Antiserum tested								
		O9	O12	O poly A-S	Anti-serogroup D	C-Se	C-St	C-Spg	P-St	C-SPF
Hb	3403	72 (4)	268 (1)	160 (2)	187 (3)	197 (6)	55 (7)	1753 (5)	76 (13)	87
Mb	4228	47 (4)	242 (3)	133 (5)	158 (4)	167 (9)	40 (12)	1674 (6)	80 (16)	107
BSA	4883	31 (11)	185 (2)	99 (16)	126 (5)	126 (11)	28 (30)	1277 (4)	68 (24)	74

^a level of immobilization.

Table 18. Biosensor responses in response units (RU) following immobilization and responses of reference antisera flowed over chip surfaces, which were prepared with **Sg LPS** oxidized in the presence of 50% (m/m) of the indicated protein. Values were corrected for the C-SPF responses, which are listed as well. Standard deviations are indicated in brackets ($N = 5$, except for the chicken and swine sera $N = 4$).

protein	Immobilization level ^a	Antiserum tested					
		O6,7	O8	O poly A-S	Anti-serogroup C	P-SI	C-SPF
Hb	5542	318 (1)	249 (11)	69 (3)	145 (10)	52 (17)	-1
Mb	9880	237 (1)	211 (12)	45 (7)	81 (26)	162 (34)	0
BSA	10344	196 (4)	124 (8)	19 (24)	61 (15)	110 (22)	-6

^a level of immobilization.

Table 19. Biosensor responses in response units (RU) following immobilization and responses of reference antisera flowed over chip surfaces, which were prepared with **SI LPS** oxidized in the presence of 50% (m/m) of the indicated protein. Values were corrected for the C-SPF responses, which are listed as well. Standard deviations are indicated in brackets ($N = 5$, except for chicken and swine sera $N = 4$).

protein	Immobilization level ^a	Antiserum tested					
		O6,7	O poly A-S	Anti-serogroup C	C-Si	P-SI	C-SPF
Hb	8180	125 (21)	16	167 (4)	377 (12)	46 (7)	-17
Mb	10819	97 (24)	8	115 (0)	349 (3)	96 (26)	-24
BSA	11364	24 (5)	-33	40 (5)	200 (2)	68 (20)	76

^a level of immobilization

- Table 20.** Biosensor responses in response units (RU) following immobilization and responses of reference antisera flowed over chip surfaces, which were prepared with St LPS oxidized in the presence of 15% (m/m) of the indicated protein. Values were corrected for the C-SPF responses, which are listed as well. Standard deviations are indicated in brackets ($N = 5$, except for chicken and swine sera $N = 4$).

protein	Immobilization level ^a	Antiserum tested								
		O4	O5	O12	O poly A-S	Anti-serogroup B	C-St	C-Spg	P-St	C-SPF
Hb	3355	419 (3)	494 (7)	212 (2)	107 (1)	809 (5)	398 (14)	619 (24)	272 (5)	9
Mb	2826	353 (8)	432 (4)	176 (8)	80 (9)	722 (2)	356 (20)	532 (24)	246 (14)	17
BSA	5588	388 (7)	304 (7)	158 (16)	82 (16)	542 (4)	324 (24)	324 (16)	244 (17)	4

^a level of immobilization.

- Immobilization levels and expected reactivity of agglutination sera with *S. enteritidis* LPS (batch Se2003.1), *S. goldcoast* (batch Sg2003.2), *S. livingstone* (batch Sl2003.1), *S. typhimurium* (batch St2003.1) and *S. meleagridis* (Batch Sm2003.1) revealed a correlation with relative amount of Hb added before oxidation (see for example Figure 3). It was also found that for a salmonella serovar-specific LPS type, the optimum Hb concentration was also production batch dependent.

On guidance of maximum response of the expected antigenic profile using a panel of standard and reference control sera and on guidance of low responses from avian SPF reference sera, optimum Hb concentration was determined for each type and for each batch of LPS (Table 21). In the text, oxidized,

haemoglobin-containing LPS preparations are further referred to as LPS^{ox}-Hb preparations.

It should be noted that, in any experiment, immobilization level did not correlate with serological responses but correlated with the amount Hb that was added.

Table 21. Effect of haemoglobin and LPS concentration on the final immobilization level of LPS derived from *S. enteritidis* (Se) and *S. typhimurium* (St).

LPS concentration ($\mu\text{g/ml}$)	Se LPS		St LPS	
	10% (m/m) Hb ^a	15% (m/m) Hb	10% (m/m) Hb	15% (m/m) Hb
25				9197
62.5 ^b	2194-2065	3173		
62.5 ^b	3400-3000			
62.5 ^b			8414	
62.5 ^b	2374			
125	3940-3647			
250	3370		10645	

^a Concentration of Hb relative to LPS

^b prepared on separate sensor channels.

The effect of dilution of Se LPS and St LPS before oxidation in the presence of 10% (m/m) or 15% (m/m) Hb relative to LPS, respectively, was investigated (Table 21). These results did not clearly reveal a correlation between LPS concentration and coupling level.

Table 22. Relationship between immobilization level of LPS^{ox}-Hb and serologic responses of O antigen antiserum and avian control sera. Fields for which a serological response was expected are shaded in a green color and values within the field are underlined. Control sera were diluted in HBS-EP containing 0.5 M sodium chloride and 0.5% (m/v) carboxymethylated dextran.

LPS batch code	immobilization level (RU)	control sera (1:20, v/v)					
		O poly A-S	C-SPF	C-Se	C-St	C-Spg	C-Si
Se2003.1	2083	<u>253.3</u>	24.9	<u>1705.9</u>	122.5	<u>4950.4</u>	58.4
	716	<u>95.8</u>	14.9	<u>681.6</u>	56.4	<u>1647.4</u>	29.1
	450	<u>48.9</u>	9.3	<u>351.5</u>	31.2	<u>733.0</u>	17.0
	292	<u>17.4</u>	10.0	<u>89.3</u>	15.2	<u>195.3</u>	10.8
Sg2003.2	8709	<u>180.0</u>	27.2	151.1	53.9	136.8	<u>184.4</u>
	8017	<u>181.3</u>	29.4	149.4	50.9	137.6	<u>186.1</u>
	6075	<u>205.0</u>	29.1	193.5	48.3	129.7	<u>206.4</u>
	3793	<u>155.4</u>	25.7	144.5	39.1	98.9	<u>157.9</u>
Sl2003.1	11526	<u>83.2</u>	51.1	69.8	76.1	115.6	<u>3708.2</u>
	11329	<u>90.3</u>	54.3	74.7	76.6	121.0	<u>3832.6</u>
	8357	<u>97.7</u>	54.8	77.3	66.8	111.9	<u>4000.1</u>
	6067	<u>88.5</u>	51.7	73.0	61.3	104.4	<u>3815.1</u>
St2003.1	2831	<u>175.0</u>	32.2	466.8	<u>607.3</u>	<u>1762.3</u>	147.4
	543	<u>48.1</u>	16.9	119.9	<u>174.9</u>	<u>369.3</u>	58.1
	321	<u>26.2</u>	11.2	64.8	<u>97.5</u>	<u>179.1</u>	36.5
	217	<u>12.2</u>	12.1	18.9	<u>28.8</u>	<u>46.0</u>	15.8

LPS stability and robustness of immobilization

Reproducibility and repeatability of the immobilization of LPS and the specific response of reference sera were tested. The Se, Sg, Sl and St LPS preparations were oxidized in the presence of their corresponding optimal Hb concentration, 5 desalted, and then stored in solution at 4°C in the dark. Under identical conditions, but accounting with the variability generally observed for immobilization of molecules at a biosensor surface, immobilization levels of oxidized LPS were either comparable in the cases of the Sl (5% RSD) and St (8% RSD) batches or tended to increase in the cases of the Se (13% RSD) and 10 Sg (8% RSD) batches over two months of storage (Table 23 through Table 26). The responses of reference sera were probed on the prepared biosensor chips as well. Inspection of these did not reveal a correlation between immobilization level and specific response. For example, immobilization levels of Se LPS and Sg LPS may tend to increase over time; this increase was not 15 reflected in the response of the binding of serum antibodies to the immobilized antigens. The relative standard deviation varies between 12% and 38%. When considering the first 3 measuring days (day 0, day 7 or 10, and day 14 or 17) the variance is greatly reduced from 3.5% to 23% (results not shown). Tables 27 to 30 summarizes the repeatability of the method at several 20 moments over a 12-months period. At each analysis time point, a fresh aliquot of Hb-fortified LPS was oxidized, immobilized and analysed and thus reflects the sum of variability of several steps.

Table 23. Response of agglutination and reference anti-sera with LPS^{ox}-Hb Se2003.1 prepared at day 0 and stored at 5-8°C. The LPS preparation was immobilized and tested after oxidation at the days indicated.

Day of analysis	immobilization level (RU)	Se2003.1 oxidized in presence of 15% Hb				
		O9	O12	O poly A-S	C-Se	C-Spg (1:100)
0	2708	446	302	330	2597	3535
10	3007	352	258	252	2205	3726
17	3194	472	320	292	2520	3586
31	3642	509	356	340	1486	3727
62	3640	361	223	156	1676	2329
Average	3238	428	292	274	2097	3381
St. dev.	407	69	52	74	498	594
RSD (%)	13	16	18	27	24	18

Table 24. Response of agglutination and reference anti-sera with LPS^{ox}-Hb Sg2003.2 prepared at day 0 and stored at 5-8°C. The LPS preparation was immobilized and tested after oxidation at the days indicated.

Day of analysis	immobilization level (RU)	Sg2003.2 oxidized in presence of 50% Hb		
		O6,7	O8	O poly A-S
0	7262	668	465	117
7	9972	714	531	110
14	10597	657	566	111
20	11453	857	502	100
28	11866	1085	504	130
59	12002	484	364	58
Average	10525	744	489	104
St. dev.	869	226	77	27
RSD (%)	8	30	16	26

Table 25. Response of agglutination and reference anti-sera with LPS^{ox}-Hb Sl2003.1 prepared at day 0 and stored at 5-8°C. The LPS preparation was immobilized and tested after oxidation at the days indicated.

Day of analysis	immobilization level (RU)	Sl2003.1 oxidized in presence of 50% Hb	
		O6,7	O poly A-S
0	11509	318	84
7	13581	289	63
14	14870	250	54
20	13582	475	65
28	14976	468	75
59	14744	202	31
Average	13877	334	62
St. dev.	707	127	17
RSD (%)	5	38	27

Table 26. Response of agglutination and reference anti-sera with LPS^{ox}-Hb St2003.1 prepared at day 0 and stored at 5-8°C. The LPS preparation was immobilized and tested after oxidation at the days indicated.

Day of analysis	immobilization level (RU)	St2003.1 oxidized in presence of 15% Hb						
		O4	O5 (1:200)	O12	O poly A-S	C-St	C-Se	C-Spg (1:100)
0	5391	594	540	278	362	441	544	1634
10	4432	463	567	242	273	448	356	1442
17	4473	517	578	250	262	362	396	1294
31	4769	559	582	296	326	440	255	1464
62	4996	442	321	209	196	340	297	983
Average	4812	515	518	255	284	406	370	1363
St. dev.	397	64	111	34	64	51	112	244
RSD (%)	8	12	22	13	22	13	30	18

Table 27. Responses of freshly oxidized **Se LPS** (batch Se2003.1) on indicated time points (in months). The LPS was isolated from bacterial cells, fortified with 15% (m/m) Hb, dried and stored at 4-7°C until day of oxidation, immobilization and analysis.

Analysis (month)	immobilization level (RU)	Se2003.1 oxidized in presence of 15% Hb				
		O9	O12	O poly A-S	C-Se (1:200, v/v)	C-Spg (1:100, v/v)
0	2708	446	302	330	2597 ^a	3535 ^b
1 ^c	1326	252	165	172	464	1274
2	2380	372	262	230	502	1478
3	2003	398	208	230	676	1864
5	2309	144	207	337	534	1539
7	3487	240	444	547	786	2368
9	3721	215	407	576	705	2415
12	1450	78	145	219	159	1272

5 ^a serum was diluted 1:50 (v/v)

^b serum was diluted 1:100 (v/v)

^c LPS was diluted at another volume ratio

Table 28. Responses of freshly oxidized **Sg LPS** (batch Sg2003.2) on indicated time points (in months). The LPS was isolated from bacterial cells, fortified with 50% (m/m) Hb, dried and stored at 4-7°C until day of oxidation, immobilization and analysis.

Analysis (month)	immobilization level (RU)	Sg2003.2 oxidized in presence of 50% Hb		
		O6,7	O8	O poly A-S
0	7262	668	465	117
1	7428	929	451	100
2	9023	639	459	91
3	13152	474	606	119
5	8087	549	446	294
7	9724	622	382	286
9	8870	692	508	364
12	7088	491	.*	281

5 **Table 29.** Responses of freshly oxidized **Sl LPS** (batch Sl2003.1) on indicated time points (in months). The LPS was isolated from bacterial cells, fortified with 50% (m/m) Hb, dried and stored at 4-7°C until day of oxidation, immobilization and analysis.

Analysis (month)	immobilization level (RU)	Sl2003.1	
		O6,7	O poly A-S
0	11509	318	84
1	11442	417	62
2	12280	287	55
3	13152	231	65
5	11896	217	159
7	11563	271	149
9	10882	298	191
12	10138	204	131

Table 30. Responses of freshly oxidized St LPS (batch St2003.1) on indicated time points (in months). The LPS was isolated from bacterial cells, fortified with 15% (m/m) Hb, dried and stored at 4-7°C until day of oxidation, immobilization and analysis.

Analysis (month)	immobilization level (RU)	St2003.1 oxidized in presence of 15% Hb					
		O4	O5 (1:200, v/v)	O12	O poly A-S	C-St	C-Spg (1:200, v/v)
0	5391	594	540*	279	362	441	1634 ^a
1 ^b	3487*	339	470	181	207	352	597
2	4623	482	325	257	242	411	612
3	4079	425	608	187	235	982	624
5	6873	359	369	167	302	393	605
7	5410	681	734	334	453	559	837
9	4008	638	664	278	441	524	824

5 ^a serum was diluted 1:100 (v/v);

^b Following oxidation, LPS-containing solution was diluted twice instead of once.

Example 2

Introduction

5 *SPR biosensor for detection of egg yolk antibodies reflecting Salmonella enteritidis infections*[AAB4]

Salmonella is one of the major causes of bacterial gastro-enteritis of humans (Fischer, 2004; van Duynhoven *et al.*, 2005). In the Netherlands, between 1994-1998, *Salmonella enterica serovar enteritidis* (*S.e.*) was the most often
10 isolated serovar (Pelt *et al.* 1999). Within this serovar, eggs and egg products were the most important source of infection. Despite several control measures, approximately 9 % of the Dutch layer flocks become infected annually. As egg contamination with *Salmonella* continues to be a threat for public health, it is important to detect an infection of a flock as soon as possible by an adequate
15 surveillance programme.

The current Dutch monitoring system in layer finisher hens is based on serology (Bokkers, 2002). The aim is to reduce the prevalence of *S.e.* and *S. typhimurium* in the layer sector. Sampling, however, occurs only twice: before and at the end of the laying period. The current surveillance programme,
20 therefore, cannot detect all infections of flocks during the layer period, and farmers cannot 'guarantee' that their products are from *Salmonella*-free layers.

Consequently, the surveillance programme should be improved. As an alternative to current serology, testing of eggs for antibodies could be
25 performed. Egg sampling has the advantage that it can be performed on egg packing plants, in a high sampling frequency and with large sample sizes.

Tests for detection of antibody in eggs have been developed and used before. The existing tests are often based on enzyme-linked immunosorbent assays (ELISA) using different (combinations of) antigenic components of *Salmonella*
30 spp. (see for examples Refs. Gast *et al.*, 2002 and 1997; Skov *et al.*, 2002; Holt *et*

al, 2000; Desmidt *et al.*, 1996; Sachsenweger *et al.*, 1994; Van Zijderveld *et al.*, 1992). Recently, the possible suitability of biosensors for the detection of humoral response has been recognized (Bergwerff and van Knapen, 2006, accepted for publication; Bergwerff and van Knapen, 2003; Jongerius-Gortemaker, 2002; Pyrohova *et al.*, 2002; Vetcha *et al.*, 2002; Li *et al.*, 2002; Liu *et al.*, 2001; Uttenthaler *et al.* 1998). A biosensor consists of a re-usable immobilized biological ligand that 'senses' the analyte, and a physical transducer, which translates this phenomenon into an electronic signal (Jongerius-Gortemaker *et al.*, 2002). The use of biosensors promises the possibility of high throughput analyses, and also the detection of multiple serovars or serogroups within a family of infectious disease agents -or antibodies against these agents- in a single run. This offers the opportunity to improve surveillance programmes, as more samples can be tested in a higher frequency during the layer period.

15 This example evaluates the sensitivity, specificity and discriminatory capacity of a surface plasmon resonance (SPR) biosensor (biacore 3000) antibody detection test in egg yolk based on the lipopolysaccharide (LPS) of *Salmonella enterica serovar enteritidis* and compares the results to those obtained with a *g,m flagellin* based commercial ELISA test kit and a LPS based commercial ELISA test kit for detection of egg-antibodies by creating and analyzing receiver operating characteristic (ROC) curves.

20

2. Materials and methods

2.1 SPR biosensor method

25 We adapted the surface plasmon resonance (biacore 3000 by Biacore AB, Uppsala, Sweden) detection method of serum antibodies against *Salmonella enteritidis* using LPS antigen as developed by Bergwerff *et al.* (in prep) to egg yolk. The first adjustment to this method was made in the sample preparation. After separation of egg yolk and egg white, the egg yolk was diluted 1:5 (v/v) in 10 mM HEPES buffer at pH 7.4, containing 3 mM EDTA, 0.15 M sodium

30

hydrochloride, 0.005 % (v/v) surfactant P20 (Biacore AB, Sweden), and additional 0.85 M sodium chloride (Merck, Darmstadt, Germany), 1 % (m/v) carboxymethylated dextran (Fluka Chemie, Buchs, Germany) and 0.05 % (v/v) Tween 80 (Merck, Germany). It was mixed with glass pearls, centrifuged at
5 15,000 *g* at ambient temperature for 25 min; the supernatant was filtrated over a 0.45- μ m filter (Schleicher & Schuell, Dassel, Germany).

A second adaptation was the cleaning of the sensor chip. Following analysis of each series of 15 egg yolk samples, a solvent containing 0.5 % (w/v) sodium dodecyl sulphate (Biacore AB, Sweden) was injected to remove deposited egg
10 yolk components.

2.2 Reference sera and egg yolks

Sera were used as reference in the various tests due to unavailability of sample stock of well-defined reference egg yolks. Lyophilized, defined SPF and
15 reference sera originating from chickens infected with a) *S. enteritidis*, b) *S. typhimurium*, c) *S. infantis*, or d) *S. pullorum* were obtained from the Animal Health Service Ltd. (Deventer, Netherlands). These sera were prepared from pooled sera. Before use, lyophilized sera were reconstituted in 1 ml Milli-Q. Additionally, monoclonal mouse anti-*Salmonella* antibody anti-group B, -group
20 C, -group D and -group E was used (Sifin, Berlin, Germany).

Internal-control-egg yolks were used to establish analytical sensitivity and repeatability of the biosensor assay. They consisted of a specific pathogen free (SPF) egg yolk sample (Animal Health Service Ltd., Netherlands) and a highly immuno-responsive pre-ovulatory follicle sample originating from experiment
25 2 (*cf.* section 2.4. below).

2.3 ELISA

Samples were assayed using sandwich enzyme immunoassay techniques. Two commercially available *S.e.* antibody detection kits were used; Flockscreen

S.e. Guildhay (Guildford, England) and FlockChek *S.e.* IDEXX (Westbrook, Maine, USA). The samples were analyzed according to the company's procedures.

The Guildhay *S.e.* indirect ELISA is based on LPS as antigen. The wells of microtiter plates were coated with LPS, 1:500 dilutions of samples were added in mono. Test results were expressed as an S / P ratio according to the following formula:

$$S / P = \frac{(\text{optical density sample} - \text{optical density negative controls})}{(\text{optical density positive controls} - \text{optical density negative controls})} \quad (1)$$

The S / P ratio was interpreted using the following criteria: Egg yolk: S / P ≤ 0.08 = immuno-negative; 0.08 < S / P < 0.25 = immuno-suspect; S / P ≥ 0.25 = immuno-positive.

The IDEXX *S.e.* competitive ELISA is based on *g,m flagellar* antigen. The wells of microtiter plates were coated with *g,m flagellar* antigen, 1:2 dilutions of samples were added in mono. The results were expressed as S / N ratio as follows:

$$S / N = \frac{\text{optical density sample}}{\text{optical density negative controls}} \quad (2)$$

The S / N ratio was interpreted using the following criteria: Egg yolk: S / N ≥ 0.75 = immuno-negative; 0.75 < S / N < 0.59 = immuno-suspect; S / N ≤ 0.59 = immuno-positive.

2.4 Experiments

The egg samples used in this study originated from two infection experiments.

Experiment 1. Fifteen one-week-old layer hens (Isa Brown) were housed in negative pressure high-efficiency particulate air filter (HEPA) isolators with a volume of 1.3 m³ and fitted with a wire floor of 1.1 m², and applying a 12 h light to 12 h dark photoperiod rhythm. The isolators were ventilated at a rate of approximately 30 m³/h. During the growing period, no *Salmonella* could be cultured from bedding. The chickens were provided with non-medicated feed and water *ad libitum*. They were housed, handled and treated following approval by the institutional animal experimental committee of the Dutch Animal Health Service Ltd. in accordance with the Dutch regulations on experimental animals. All hens were inoculated orally once with 1 x 10⁸ CFU per bird in week 20 of the experiment using *S. enteritidis* CL344 (Animal Health Service Ltd., Deventer, The Netherlands). Before and after inoculation, eggs were collected on a daily basis, but not labeled individually and not dated. The eggs were stored at ambient temperature for four weeks and subsequently at 4 °C. The experiment ended in week 22 and produced 147 'positive' and 71 'negative' samples.

Experiment 2. This experiment is described in detail by Van Eerden *et al.* (2005, in prep). In short, 128 15-week-old layer hens (Lohmann Brown, 16 birds, 8 replications) were divided into two groups (8 hens each) and housed individually in two climate cells, used as isolators, in the same room, under a 9 h light and 15 h dark photoperiod rhythm. During the growing period until inoculation, no *Salmonella* was cultured from feces. They were provided with non-medicated feed and water *ad libitum*. The animal experiment was conducted according to the Guidelines for Animal Experimentation of Wageningen University and approved by the Ethical Committee under Reference Number of 2003219. Each of sixty-four hens (16 birds, 4 replications) was inoculated orally once with 1 x 10⁸ CFU nalidixic acid-resistant *S.e.* (ASG, Lelystad, the Netherlands) one week after the experiment started. The other 64 birds were considered uninfected controls. Eggs were collected at day 21 and day 28 post-inoculation. Twelve times, the eggs from

one climate cell were collected and pooled after cracking of the shell (cf. section 2.5 below). Four of the pooled egg samples were taken from 'uninfected' climate cells. Besides the collection of pooled egg samples, ten egg samples were taken from ten individual birds, of which seven were uninfected. The egg yolk and white were mixed and stored at -20°C . The experiment ended four weeks after inoculation. Pre-ovulatory follicles were then harvested from eight uninfected and five infected individual birds.

2.5 Preparation of egg samples

Eggs from experiment 1 were prepared in the following manner. To facilitate aseptic preparation, the eggshells were disinfected with a 70 % (v/v) aqueous ethanol. Subsequently, the eggs were cracked, and the contents were collected in sterile petri dishes. A volume of 1 ml egg yolk was collected using a sterile disposable syringe and portioned in 200 μl fractions. Each fraction was diluted with a buffer appropriate for either *SPR* biosensor or ELISA analysis, and then stored at -20°C .

Likewise, pooled egg yolk and white and pre-ovulatory follicles obtained from experiment 2 were fractionated, diluted and stored at -20°C .

2.6 Evaluation of the SPR biosensor method

2.6.1 Analytical sensitivity and specificity

To establish the limit of detection of the assay, eight 1:2 (v/v) serial dilutions of a highly immuno-responsive egg yolk sample and a SPF negative control sample were analyzed in triplicate by the *SPR* biosensor. Analysis of variance (ANOVA) was performed using SPSS (SPSS for Windows, Standard Version, 1999) to evaluate differences in *SPR* biosensor responses obtained after injection of the serially diluted control samples.

Reference sera were used to spike a SPF yolk for to test the specificity of the *SPR* biosensor assay. For this purpose, egg yolk was spiked with 1) *S*.

enteritidis- (serogroups D), 2) *S. infantis*- (serogroups C), 3) *S. pullorum*- (serogroups D), and 4) *S. typhimurium*- (serogroups B) reacting antisera.

These samples were diluted by their volumes either at a rate of 1:100 (1, 2 and 3) or at 1:50 (4). Further specificity testing was performed by spiking SPF egg
5 yolk with 1:100 (v/v) diluted mouse monoclonal antibody reacting with *Salmonella* serogroups B, C, D and E.

2.6.2 Repeatability

The repeatability of the *SPR* biosensor assay was assessed by running the
10 highly immuno-responsive egg yolk sample and the SPF negative control egg yolk sample twice on a single day and on three consecutive days (in triplo). Means, standard deviations (SD) and percent coefficient of variation (%CV) values were calculated in Excel 2000 (Microsoft software package).

15 2.6.3 ROC curves

Receiver operator characteristic (ROC) curves were generated using the results from the *SPR* biosensor and ELISA analyses to assess the test performances of each assay (Zweig and Campbell, 1993). Using SPSS, the overall accuracy of each assay was calculated from the integrated area under
20 the curve (AUC), corresponding standard error (SE) and the probability of the null hypothesis of the true AUC being 0.5. By use of non-parametric ROC analysis (Metz *et al.*, 1998), the accuracy of *SPR* biosensor assay detection of antibodies against *S.e.* was compared with the accuracy of the two ELISA's. The gold standard was the infection status of the experimental group.

25

2.6.4 Diagnostic sensitivity and specificity

In a ROC curve the true positive rate (sensitivity) is plotted in function of the false positive rate (100-specificity) for different cut-off points of a parameter. Each point on the ROC curve represents a sensitivity/specificity pair
30 corresponding to a particular decision threshold. Thus, the maximum

diagnostic sensitivity at the highest diagnostic specificity for the *SPR* biosensor assay and the two ELISA's were calculated, using SPSS. For the *SPR* biosensor test using the samples from experiment 1, the maximum diagnostic specificity at the highest diagnostic sensitivity and the optimal
5 combined diagnostic sensitivity and specificity were also calculated.

3. Results

3.1 Analytical sensitivity and specificity

10 A 1:640 (v/v) dilution of the highly immuno-responsive egg yolk sample was, at 50 RU, the highest dilution tested that differed significantly ($P < 0.001$) from the negative control.

The test signal of the SPF egg yolks spiked with *S. enteritidis*- (1:100, 145 RU), *S. pullorum*- (1:100, 1012 RU) or *S. typhimurium*- (1:50, 58 RU) positive
15 sera were above the optimized cut-off value of 52 RU (cf. section 3.4.1 below) and considered positive, as was the SPF egg yolk spiked with mouse anti-*Salmonella* group D (1:100, 130 RU). Non-spiked SPF yolk was found to be negative, i.e. average response was 30 RU. The yolks spiked with *S. infantis*- (1:100, 24 RU) positive serum and mouse antiserum against *Salmonella*
20 serogroups B (1:100, 27 RU), C (1:100, 16 RU), and E (1:100, 15 RU) were also below the cut-off value.

3.2 Repeatability

The coefficient of variation within a single day was 1 % for the highly
25 immuno-responsive egg yolk sample and 13 % for the negative sample. The coefficient of variation from day-to-day during three days was 2 % for the positive sample and 17 % for the negative sample.

3.3 Threshold determination

30 3.3.1 ROC analysis

ROC analysis was performed on the assay results of 71 and 135 egg yolk samples from uninfected and infected chickens, respectively, from experiment 1 (not all tests were performed on 12 samples from infected chickens).

Integrated areas under ROC curves were 0.892 (SE 0.024, $P < 0.001$) for the SPR biosensor assay; 0.432 (SE 0.039, $P = 0.103$) for the IDEXX ELISA and 0.430 (SE 0.039, $P = 0.096$) for the Guildhay ELISA (Table 31). The ROC curves are depicted in Fig. 4. The integrated area (AUC), and thus the overall accuracy, for the SPR biosensor assay was significantly larger than those of the IDEXX ($Z = 11.5$, $P < 0.001$) and Guildhay ELISA ($Z = 10.5$, $P < 0.001$).

ROC analysis was also performed for four combined egg white and yolk samples and 15 egg yolk samples from uninfected, and eight combined egg white and yolk samples and eight egg yolk samples infected chickens from experiment 2. The integrated areas under ROC curves were 0.811 (SE 0.082, $P = 0.002$) for the SPR biosensor assay, 0.615 (SE 0.098, $P = 0.098$) for the IDEXX ELISA and 0.870 (SE 0.064, $P < 0.001$) for the Guildhay ELISA (Table 32 and Fig. 5). The AUC of the SPR biosensor assay was significantly ($Z = 1.9$, $P = 0.055$) larger than that of the IDEXX ELISA, but not different from that of the Guildhay ELISA ($Z = -1.0$, $P = 0.322$).

3.4 Performance estimates

3.4.1 Diagnostic sensitivity and specificity estimates

With respect to the results of the samples acquired from Experiment 1, samples from the uninfected population gave biosensor responses ranging from 6 to 50 RU. The responses of the samples from the infected population ranged from 11 to 3584 RU. At a cut-off value of 52 RU, 24 out of 135 samples had to be considered immuno-negative.

A cut-off value of 52 RU yielded the highest possible diagnostic specificity estimate of 100 % (with a 95 % exact confidence interval (CI) of 95-100 %) and a diagnostic sensitivity estimate of 82 % (95 % CI: 76-98 %) for the SPR

biosensor assay test. A cut off value of 10 RU yielded the highest possible

diagnostic sensitivity estimate of 100 % (95 % exact CI: 97-100 %) and a - specificity estimate of 1 % (95 % CI: 0-4 %). A cut-off value of 42 RU yielded the optimal combined diagnostic sensitivity and -specificity: 84 % (95 % CI: 77-90 %) and 99 % (95 % CI: 96-100 %), respectively.

- 5 At a cut-off value of OD_{550nm} 0.11, the IDEXX ELISA had a diagnostic specificity of 100 % and a -sensitivity of 1 % (95 % CI: 0-3 %). The OD_{550nm} of the samples from the uninfected population ranged from 0.174 to 1.377, i.e. in excess of the cut off value at 0.11. Of the positive population, 145 out of 147 samples had to considered immuno-negative at the chosen cut-off value,
- 10 namely corresponding OD_{550nm} ranged from 0.042 to 1.572. The Guildhay ELISA had a diagnostic specificity of 100 % and a -sensitivity of 16 % (95 % CI: 10-22 %) at a cut-off value of OD_{650nm} 0.12. None of the samples from the uninfected population showed OD_{650nm} values in excess of 0.12 (0.051 to 0.093). In case of the positive population, 124 out of 147 samples had to be considered
- 15 immuno-negative. The OD_{650nm} of these samples ranged from 0.048 to 1.471.

- In the case of Experiment 2, a cut-off value of 542 RU yielded the highest possible diagnostic specificity estimate of 100 % (95 % exact CI: 82-100 %) and a diagnostic sensitivity estimate of 63 % (95 % CI: 39-86 %) for the *SPR* biosensor assay test. The samples from the uninfected population had RU
- 20 values ranging from 101-448. The infected population values ranged from 117-3012 and 6 out of 16 samples had negative test results. At a cut-off value of OD_{550nm} 0.49, the IDEXX ELISA had a diagnostic specificity of 100 % and a -sensitivity of 19 % (95 % CI: 0-38 %). None of the samples from the uninfected population had OD_{550nm} values of less than 0.49. The values ranged from 0.537-
- 25 1.621. Of the positive population, 13 out of 16 samples had negative test results at the chosen cut-off value. The values ranged from 0.134-1.630. The Guildhay ELISA had a diagnostic specificity of 100 % and a -sensitivity of 67 % (95 % CI: 43-91 %) at a cut-off value of OD_{650nm} 0.14. None of the samples from the uninfected population had OD_{650nm} values of more than 0.14. The values
- 30 ranged from 0.072-0.140. Of the positive population, 5 out of 15 samples had

negative test results (one sample could not be tested). The values ranged from 0.086-2.144.

4. Discussion

5 The aim of this study was to quantify the test characteristics of the *SPR* biosensor for the detection of *S.e.* antibodies in eggs. The results showed that the *SPR* biosensor assay performed significantly better than the two commercially available ELISA's for samples from Experiment 1. The combined optimal diagnostic sensitivity and -specificity of the *SPR* biosensor was 84 %
10 (77-90 %) and 99 % (96-100 %), respectively. Neither the *g,m flagellin*-based IDEXX ELISA, nor the LPS-based Guildhay ELISA were able to detect *S.e.* infection with a higher combined diagnostic sensitivity and specificity using this test panel. This study indicates that an *SPR* biosensor assay could be a new and powerful tool for monitoring *Salmonella enterica serovar enteritidis*
15 infections in layer flocks through antibody detection in eggs.

The *SPR* biosensor assay offers the possibility of detecting infections in a fast and reliable way. The high quality of the test and the technical and animal welfare advantages of egg collection are good reasons to explore its use for screening of populations. In addition, the configuration of the applied *SPR*
20 biosensor from Biacore allows the simultaneous detection of antibodies to multiple *Salmonella* serovars in a single run in a single sensor channel or in separate sensor channels on the same sensor chip (results not shown). This could be of significance because it is well known that serovars differ over countries and over time (see for examples Refs. Guerin *et al.*, 2005; van
25 Duijnkeren *et al.*, 2002).

The test evaluation was carried out using eggs from two experiments that were not carried out specifically for this test evaluation, possibly influencing test performance. The 'positive' eggs were collected probably at a time point that humoral response was developing in the exposed chickens. These
30 'premature' eggs were analyzed and their false-negative results interfere with

the evaluation of the assays. Could it have been possible to exclude eggs until 2 weeks post-infection, the diagnostic sensitivity of each test, ELISA or SPR biosensor, would have been improved.

Antibody detection in serum is more sensitive than in eggs, because the
5 appearance of antibodies in eggs is preceded by the appearance in serum by a week (Gast and Beard, 1991; Sunwoo *et al.*, 1996; Skov *et al.*, 2002). However, flock sensitivity of tests for antibodies in eggs can be improved by taking more samples, which is easier when using eggs.

The biosensor performance (AUC 0.892) was compared to that of two
10 commercial ELISA's, (IDEXX AUC 0.432, Guildhay AUC 0.430). To our knowledge the IDEXX ELISA was not validated for eggs, but quantitative data exist about the test's performance in comparison to other tests: Van Zijderveld *et al.* (1992) evaluated four different ELISA's for diagnosis of *S.e.* infections in experimentally infected chickens. They reported a specificity of 100 % and a
15 sensitivity of 95 % for 127 egg yolks from eggs laid between 13 and 40 days after infection with *S.e.*. In our evaluation, the IDEXX test performed not as well as in the 1992 evaluation. An explanation could be the different sample selection, as our samples originated from infection experiments that stopped at 2 and 4 weeks after inoculation, having had less time to develop a humoral
20 response.

Shared O-antigens among members of *Salmonella* serogroups B and D are known to limit the specificity of detecting *S.e.* using lipopolysaccharide antigens (de Vries *et al.*, 1998; Baay and Huis in 't Veld, 1993; Hassan *et al.*, 1990). This is confirmed by our results: the assay could not differentiate
25 between infections with serovars *enteritidis*, *gallinarum* and *typhimurium*, sharing O 9 and O 12. As the zoonotic serovars of the three (*S. typhimurium* plus *S. enteritidis*) represent 80% of isolates identified by the national reference laboratories participating in the Enter-net surveillance network between 1998-2003 (Fischer *et al.*, 2004), this finding has limited clinical
30 relevance for the human population. The assay did differentiate between SPF

egg yolk spiked with mouse anti-*Salmonella* group B (1:100, 27 RU) and D (1:100, 130 RU). This is not surprising, because the LPS of *Salmonella enteritidis* has O 1, O 9 and O 12 as somatic antigens, whilst the group specific test reagents contain the following monoclonal antibodies; anti-*Salmonella* group B: Anti-O 4, O 5, O27; anti-*Salmonella* group D: Anti-O9.

The cut-off value from Experiment 2 was much higher than the cut-off from Experiment 1, possibly because part of our samples consisted of egg white and yolk instead of egg yolk only.

For different applications, different cut-off values may be optimal. Relative costs or undesirability of errors (false positive/false negative classifications) and the expected relative proportions of infected and uninfected hens are important parameters in the determination of the cut-off value, which affects the diagnostic value of the assay. We would suggest a cut-off value which minimizes the number of false positive results, reasoning that frequent sampling and testing would be necessary if the assay was to be used in a surveillance programme in the layer population, given the relatively low prevalence of *S.e.*.

The *SPR* biosensor technique has successfully detected egg antibodies to determine experimental infections in chickens. In future screening programmes, the *SPR* biosensor could possibly detect different analytes at the same time.

Table 31. ROC analysis of the results of samples derived from Experiment 1 analyzed by *SPR* biosensor, IDEXX and Guildhay ELISA's.

Characteristic	<i>SPR</i> biosensor assay	IDEXX ELISA	Guildhay ELISA
Optimized cut-off	52 RU	OD _{550nm} 0.11	OD _{650nm} 0.12
Diagnostic sensitivity (%)	82	1	16
95% CI (%)	76-89	0-3	10-22
Diagnostic specificity (%)	100	100	100
95 % CI ^a	95-100	95-100	95-100
AUC	0.892	0.432	0.430
95 % CI	0.844-0.939	0.356-0.508	0.355-0.506

^a Fisher's exact test

Table 32. ROC analysis of the results of samples derived from Experiment 2 analyzed by *SPR* biosensor, IDEXX and Guildhay ELISA's .

<i>Characteristic</i>	<i>SPR</i> biosensor assay	IDEXX ELISA	Guildhay ELISA
Optimized cutoff	542 RU	OD _{550nm} 0.49	OD _{650nm} 0.14
Diagnostic sensitivity (%)	63	19	67
95% CI	39-86	0-38	43-91
Diagnostic specificity (%)	100	100	100
95 % CI ^a	82-100	82-100	82-100
AUC	0.811	0.615	0.870
95 % CI	0.649-0.972	0.424-0.806	0.745-0.996

^a Fisher's exact test

Example 3

Introduction

Direct detection of Campylobacter spp. through monitoring bacteriophage infections using LPS-coated beads

- 5 *Campylobacter* is the most commonly food-borne pathogen in developed countries, causing gastroenteritis characterized by watery and/or bloody diarrhea. *Campylobacter* is associated with Guillain- Barré (GBS), Reiter's and haemolytic uremic (HUS) syndromes and reactive arthritis (FSAI, 2002; Lake et al., 2003; Tauxe, 2000). In the last 20 years, the infection rate of
- 10 *Campylobacter* is still increasing in many developed countries, maybe due to the improvements in detection and reporting. In the United States of America, 2,400,000 cases of campylobacteriosis are reported annually corresponding to approximately 1% of the USA population (Tauxe, 2000).

- Wild birds and domestic animals are reservoirs for *Campylobacter* and shed
- 15 bacteria to the environment. Poultry is an importance vehicle for *Campylobacter* infection in humans. Indeed, strains, which were isolated from chickens, could be isolated from patients as well (Coker, 2000).

- Epidemiological studies have shown that consumption and handling of poultry meat should be considered as a major risk for human infection with *C. jejuni* or
- 20 *C. coli* (FSAI, 2002). The most consistent risk factor in United States, New Zealand and Europe has been consumption or contact with raw or undercooked poultry, accounting for 10% to 50% of all cases of campylobacteriosis (Tauxe, 2000). *C. jejuni*, *C. coli* and *C. lari* represent about 90% of human campylobacteriosis (Stern and Line, 2000). The infective dose of
- 25 *Campylobacter* is considered to be low, ranging from 500 – 10,000 cells (FSAI, 2002).

Campylobacter are Gram negative, curve, S-shaped, or spiral shaped bacilli having one or two flagella at one of the poles and highly motile (Christensen et

al., 2001). *Campylobacter* grows between 30.5°C and 45°C at an optimum temperature of 42°C. Optimum growth is established at 10% carbon dioxide, 5-6% oxygen, and 85% nitrogen (FSAI, 2002).

Traditional phenotyping methods for determination of *Campylobacter* take 4 to 5 days and involve pre-enrichment followed by isolation from selective agar and confirmation by biochemical test. Due to the perishable nature of food items and the speed required for analysis of food products more rapidly, sensitive and specific methods are needed for cost-effective *Campylobacter* detection.

Immunomagnetic separation (IMS) procedures were used by Waller and Ogata (2000), Che et al. (2001), Yu et al. (2001) to concentrate *C. jejuni* from poultry meat without pre-enrichment cell culture step. This approach could retrieve 10⁴ colony forming units (cfu)/g in poultry meats as detected with atomic force and fluorescence microscopy (Yu et al., 2001). IMS can potentially reduce pre-enrichment time of *Campylobacter* and may overcome the problems of inhibitors from food sources such as PCR inhibitors (Benoit and Donahue, 2003). The use of IMS may thus speed up the enrichment of the analyte.

This example describes a down-stream detection method using *Campylobacter*-specific bacteriophages, i.e. small viral organisms that attach to or infect living *Campylobacter* bacteria. Their attachment or infection is dependent of the phase of life cycle of the bacterium. Binding to or infection of *Campylobacter* may namely occur in the stationary, log or lag phase of the bacterium and depends of the phage species as well. Infection of the bacterium results usually in a high number of copies of the bacteriophage. Recording this increment of phages is therefore used as an analytical instrument to trace the presence of *Campylobacter* in the original sample.

The aim of this study is to demonstrate the application of bacteriophages as specific and sensitive analytical tools for the detection of *Campylobacter* in animal products, such as faeces and (poultry) meats.

Materials and methods

Experimental set-up

Following homogenisation, e.g. facilitated by stomaching, IMS will be used to purify and concentrate *Campylobacter* from contaminated samples, such as meat and faeces. In a second step, IMS-isolated bacteria will be incubated with an appropriate strain of bacteriophage. Non-attaching bacteriophages will be washed from the cell isolate using the same IMS procedure. Infected and/or bacteriophage-carrying IMS-immobilised *Campylobacter* are then introduced in a fresh and pure culture of reference *Campylobacter* that is in a stationary phase. This cell culture is used as a foreign host to boost the multiplication of the bacteriophages. Following a short culture to allow the bacteria to reach their log-phase, bacteriophages will be harvested by centrifugation. The bacteriophage-containing supernatant will be incubated with LPS-coated fluorescent beads. Here, the bead is coated as described in the Example with the LPS isolated from *Campylobacter* used as the host organism. The presence of bacteriophages bound to the fluorescent beads will be tested in two ways. Following the addition of and incubation with anti-bacteriophage antibodies tagged with a fluorescent label, the amount of fluorescence will correspond with the concentration of bacteriophages and indirectly with the concentration of campylobacter in the original sample. In an alternative approach, anti-LPS antibodies containing a fluorescent tag will compete with bacteriophages for binding places. A decrease of recorded fluorescence compared to a *Campylobacter*-free sample will, therefore, indicate a *Campylobacter* positive sample.

The test will be validated in terms of selectivity and sensitivity for *C. jejuni*, *C. coli* and *C. lari* in different matrices, including faeces, skin and meat from pigs and chickens. Closely related organisms, such as *Arcobacter* species, will be used to test the specificity of the method.

Bacterial and viral strains and culture condition

C. jejuni (ATCC 33291) and *C. coli* (ATCC 33559) will be bought from Microbiologics (St. Cloud, USA). The bacteria will grow in tryptone soya broth (TSB) (Oxoid, CM 129, Hampshire, England) for 24 h at 42°C, under
5 microaerophilic atmosphere, which will be generated using a gas package (BBL, Becton Dickinson, Sparks, USA). *Campylobacter* are then plated onto Charcoal-Cefoperazone-Deoxycholate Agar (mCCDA) (*Campylobacter* blood-free selective agar base [Oxoid, CM 739] with CCDA selective supplement [Oxoid, SR155], cefoperazone 32 µg/ml and amphotericin B 10 µg/ml) and
10 incubated under microaerophilic atmosphere for 24 to 48 h at 42°C. One colony of pure *Campylobacter* is then transferred to tryptic soya agar (TSA) (Oxoid, CM131) and will be incubated under microaerophilic atmosphere for 24 to 48 h at 42°C and will then placed in a refrigerator at 4°C until use.
Campylobacter-infecting bacteriophages NTCC12669 , NTCC12670 ,
15 NTCC12671 , NTCC12672 , NTCC12673 , NTCC12674 , NTCC12675 , NTCC12676 , NTCC12677 , NTCC12678 , NTCC12679 , NTCC12680 , NTCC12681 , NTCC12682 , NTCC12683 , NTCC12684 are acquired from the National Type Culture Collection (London, United Kingdom).

Sample preparation

20 The pure *Campylobacter* culture stored at 4°C will be subcultured in TSB and incubated under microaerophilic atmosphere for 24 h at 42°C. This is the host for exponential growth of the bacteriophage.
An amount of 25 g of ground chicken fillet will be suspended in 225 ml of Preston broth (Nutrient broth No.2 [Oxoid, CM 67], 5% (v/v) lysed horse blood
25 [Oxoid, SR48], *Campylobacter* growth supplement [Oxoid, SR232] and modified Preston *Campylobacter* selective supplement [Oxoid, SR204]) contained by a stomacher bag. The Preston broth medium will be prepared according to the manufacturer's instruction. The sample-containing stomacher bag will be homogenized thoroughly for 90 s in a stomacher (Interscience, St.Nom,

France). The entire suspension will be then be incubated under microaerophilic atmosphere at 42°C for an appropriate incubation time to allow growth of *Campylobacter*.

5 *Immunomagnetic separation*

After enrichment with Preston broth, the stomacher bag containing the sample will be placed into the incubation pot of the IMS machine (Pathatrix™, Microscience, Cambridgeshire, UK). The apparatus is then operated according to the instructions of the manufacturer. Briefly, 50 µl of anti-*Campylobacter*
10 magnetic beads (Microscience) will be added to the sample, which is then recirculated 30 min at 37°C. The magnetically-immobilized beads are released, washed with 100 ml of pre-warmed buffered peptone water (peptone [Becton Dickinson] 10 mg/ml, sodium chloride [Merck, Darmstadt, Germany] 5 mg/ml, disodium hydrogen phosphate dihydrate [Merck] 4.5 mg/ml, potassium
15 dihydrogen phosphate [Merck] 1.5 mg/ml adjusted to pH 7.2) and then drawn to the magnet again. Wash solution was removed leaving a 200 µl suspension for selective growth and bacteriophage analyses.

Detection of bacteriophages

20 *Campylobacter*-carrying IMS-beads are contacted with a small volume of bacterium-specific bacteriophages. Following a short incubation to allow specific attachment of the phages to the surface of the targeted bacterium, IMS beads are washed and sampled to set a reference point in the final analysis procedure. The rest of the suspension is mixed with a suspension of fresh
25 *Campylobacter* species to host the growing bacteriophage. Following incubation at 42°C, the suspension is centrifuged and the supernatant will be supplemented with a volume of *Campylobacter* LPS-coated fluorescent beads. Multiplication of the phages is then assessed following the addition of either fluorescently labelled anti-bacteriophage antibodies or fluorescently labelled

anti-*Campylobacter* antibodies. Following an incubation of 15 min, the beads are analysed using e.g. a BioPlex device (Bio-Rad) to screen fluorescence immobilised on the beads as a result of specific binding reactions.

Example 4

Introduction

Detection of anti-Salmonella antibodies in porcine serum and meat juice from chickens using fluorescent beads

- 5 Micro-organisms include a wide variety of bacteria, moulds (fungi), parasites and viruses. Pathogenic micro-organisms have attracted much attention from the public as consumers of contaminated food and water, which resulted in family or community outbreaks. As a consequence, the media and politicians have played their part in increasing consumer awareness and new legislation
10 is in preparation or already in force.

- With respect to pathogenic micro-organisms, special attention is drawn to a number of zoonotic diseases, *i.e.* microbes transmissible from animals to human, for the following reasons: 1) most food- and waterborne diseases in human are zoonotic by nature; 2) many zoonotic agents have their
15 transmission route through the environment, and 3) both contamination of food/water and environment are also used by (bio)terrorists to acquire maximum impact in the society.

- Microbiological hazards can enter food chains at any point during pre-harvest, production, processing, transport, retailing, domestic storage or meal
20 preparation. From their introduction on feed or food, highly complex environments can occur in which the micro-organism can elude detection and inactivation. Efficient international distribution systems and rapid changes in consumer preferences can facilitate the swift penetration of pathogens through large populations, greatly shortening the reaction time available to public
25 health agencies.

Authorities and food producers are convinced that rapid fast, versatile and selective (diagnostic) assays are needed for environmental, feed and food monitoring to react adequately to contaminated links in the food chain. A

large portion of the explored monitoring techniques involved the use of affinity assay technologies, including biosensor platform.

In[AAB5] principle, detection of the presence of micro-organisms can be carried out in two ways: directly or indirectly. In the direct assay, the organism itself
5 is detected usually with the application of antibodies reacting with (sub)species- and/or strain-specific antigenic structures. This immunochemical analysis follows time-consuming sample preparation through culturing in selective growth media. In the case of parasite infections, this is not possible and direct detection involves microscopic inspection of samples. In the indirect
10 assay, the presence of the micro-organism is suggested by the detection of humoral (immunoglobulins) or cellular (*e.g.* cytokines) products of an immunological response of the infected host. In most studies, well-defined antigens are used to capture host's immunoglobulins in any body fluid (serology). The observed binding then reveals the nature of an invasive
15 infestation of a pathogen.

The advantages and disadvantages of indirect and direct pathogen detection are clear: *i*) individuals are not always immunologically responding to an infection; *i.e.* differences between low or high immune responders, *ii*) humoral responses are delayed several days or even weeks possibly leaving a recent
20 infection unnoticed, *iii*) serum antibodies can be found where the causative organism is not detected, as it has been rejected or retracted itself in certain (non-sampled) tissues, *iv*) serological investigations are very fast and offer better possibilities for high-throughput than direct detection, and *v*) serologic analysis of serum or plasma predicts the *Salmonella* infection status of a flock
25 or herd better than direct antigen analysis, *i.e.* classical selective bacterial culturing.

In fact, serology outperforms direct, and in most cases insensitive detection of tissue parasites, which can only be carried out by histochemistry or digestion techniques and microscopy. Significant differences are also apparent in

sample collection and preparation: whereas bacteria, fungi and viruses have to be cultured from matrices to facilitate their detection in enriched solutions, blood is relatively easily collected and prepared for analysis. Here, it should be noted, however, that antibodies can not only be retrieved from blood, plasma or serum, but also from muscle (meat juice), milk, colostrums, cerebrospinal fluid and eggs. In particular, sampling of eggs, meat juice and/or milk is easier and more cost-effective than the sampling of blood, plasma, serum or cerebrospinal fluid.

Diagnostic methods based on serologic analysis of antibody-containing biological materials are therefore supportive in so-called logistic slaughtering of animals. In this innovative processing approach, evidence-based and reliable decisions are made on the basis of continuous and intensive monitoring on farm level whether animals are allowed to enter a *Salmonella*-free or a *Salmonella*-contaminated processing infrastructure.

Among the components of the antigenic structure of the genus *Salmonella*, the somatic antigens are important as an instrument to trace immune response in animals upon an invasive infection of this organism. Somatic antigens are located on the polysaccharide part of lipid polysaccharide (LPS), which is a constituent of the bacterial cell wall. Detection a humoral response with carefully chosen LPS, the identity of the serogroup of the infecting *Salmonella* can be deduced.

In Denmark, Germany, Greece and The Netherlands, 39.5% of all *Salmonella*-positive pigs sampled at the abattoir were determined as *S. typhimurium*.

Dependent of country, other important isolates from pigs were *S. derby* (17.1%), *S. infantis* (8.0%), *S. panama* (5.1%), *S. ohio* (4.9%), *S. london* (4.4%), *S. livingstone* (3.1%), *S. virchow* (2.7%), *S. bredeny* (2.1%), *S. mbandaka* (1.1%), *S. Brandenburg* (1.0%), *S. goldcoast* (0.8%).

In case of chickens, 14% of the chickens were *Salmonella*-positive at flock level in 2002 in The Netherlands. The predominant serovar was in that case *S. paratyphi B var. java*.

The retail level a comparable percentage (13.4%) was found in the Netherlands. The most frequent *Salmonella* serovars isolated from broilers in 14 EU member states were *S. paratyphi B var. java* (24.7%), *S. enteritidis* (13.6%), *S. infantis* (8.0%), *S. virchow* (6.7%), *S. livingstone* (5.7%), *S. mbandaka* (5.5%), *S. typhimurium* (5.3%), *S. senftenberg* (5.0%), *S. hadar* (3.7%). *S. paratyphi B var. java* is dominating, but this is fully attributable to The Netherlands.

In food-producing chickens and swine, the prevalently occurring *Salmonella* serogroups are thus belonging to groups B, C and D, and in the case of swine also E.

10

In this study, a new analytical affinity assay platform is explored for the indirect detection of *Salmonella* infection in pigs and chickens. This technology platform from Luminex analyses internally coded beads which can be coated with different antigens in a single test. Only when both fluorescence of bead and bound analyte pass the detector a response will be recorded. This approach is applied to detect anti-*Salmonella* antibodies in serum and meat drip.

Materials and methods

20 Chemicals

Amine coupling kit, consisting of *N*-hydroxysuccinimide (NHS), 1-ethyl-3-(3-dimethylaminopropyl)carbodiimide hydrochloride (EDC) and ethanolamine hydrochloride – sodium hydroxide pH 8.5 were bought from Biacore AB (Uppsala, Sweden). Ethanol and trichloroacetic acid (TCA) were purchased from Merck (Darmstadt, Germany). Sodium cyanoborohydride and carbonyldiimidazole were obtained from Fluka Chemie GmbH (Buchs, Switzerland). Porcine hemoglobin (Hb) was acquired from Sigma Chemical Company (St. Louis, MO, U.S.A.). Water was obtained from of a Milli Q water purification system (Millipore, Bedford, MA, U.S.A.).

Materials

NAP-5 columns (0.5 ml; Sephadex G-25) were purchased from Amersham Biosciences and were used as described by the producer. CM5 biosensor chips were bought from Biacore AB. Dialysis bag (Spectra/Por) with a cut-off of 5 1 kDa was obtained from Spectrum Laboratories Inc. (Rancho Dominguez, CA, U.S.A.). Alexa532 was from Molecular Probes (Leiden, The Netherlands). Goat anti-swine IgG (H+L) was ordered from Jackson ImmunoResearch (West Grove, PA, USA). This antibody was conjugated with Alexa532 using standard labelling procedures.

10 *Anti-Salmonella antisera*

Salmonella monovalent 'O' somatic monoclonal antisera against O4, O5, O6₁, O7, O8, O9, O10 were purchased from Sifin (Berlin, Germany). Antibody solutions were diluted in 50 mM PBS to their working concentrations.

Reference avian and porcine sera

- 15 All reference sera were obtained from the Dutch Animal Health Service (Deventer, The Netherlands). The obtained avian reference sera were reactive with *Salmonella enteritidis* (serogroup D₁), *S. typhimurium* (serogroup B), *S. pullorum/gallinarum* (Spg; serogroup D₁) and *S. infantis*-(serogroup C₁), and were further referred to as CH-Se, CH-St, CH-Spg and CH-Si,
- 20 respectively. These chicken sera were originally prepared for ELISA analyses as positive references. In addition, specific pathogen-free chicken serum (further referred to as CH-SPF) was purchased from this organisation as a negative control reference sample. These sera were reconstitution from freeze-dried material by addition of water at a volume indicated by the producer.
- 25 Likewise, *porcine* sera from animals challenged with *S. typhimurium* and *S. livingstone* (serogroup C₁) were acquired from GD and were referenced as P-St and P-Sl, respectively. In addition, *Actinobacillus pleuropneumoniae* serotype 2-reacting *porcine* serum (GD) used as control in a Complement

Fixation Test, was exploited as negative control for porcine serum in the salmonella biosensor assay.

Methods

5 *Extraction of LPS*^[AAB6]

Overnight cultures of salmonella were prepared by applying 100 µl from their corresponding stocks on each of the 120 plates containing brain heart infusion agar (BH1a, Oxoid). The bacteria were harvested from the surface of the plates into 1 ml 9 g/l NaCl (saline) solution per agar plate using a trigalski
10 spatula. Each plate was washed twice with 2 ml saline solution. Combined bacteria were centrifuged in six tubes at 10,000 g and 4°C for 15 min and supernatant was discarded. This centrifugation step was repeated twice with 75 ml saline wash solution per tube each run. While kept on ice, pelleted bacteria were suspended in water at a volume ratio, which was a 5-fold to the
15 weight of the bacteria. An equivolume of 0.250 M (Se) or 0.500 M (Sg, Sl, Sm and St) TCA was added to give end concentrations of 0.12 M and 0.25 M, respectively, followed by continuous stirring at 4°C for 3 h. A lipopolysaccharide (LPS)-containing supernatant was then acquired at 20,000 g and 4°C for 30 min. The pH of the supernatant was adjusted to pH
20 6.5 with 5 M sodium hydroxide and when nearing the aimed pH with 0.10 M sodium hydroxide. The final volume of the LPS-containing solution was determined prior to storage at -18°C for 30 min. The solution was diluted with a double volume of freezing cold absolute ethanol from a -18°C storage place, and incubation was continued overnight at -4°C without stirring in a closed, in
25 house-build device with circulating cold ethylene glycol/water (1:4, v/v). An LPS-containing pellet was obtained after centrifugation at 20,000 g and -4°C for 30 min. The particulate material was suspended in a volume of 0.5 ml water per gram original bacterial mass weighed at the start of extraction process. The suspension was dialyzed in a 1-kDa dialysis bag against water at

4°C for two days with regular intermittent refreshment of the water. The bag content was centrifuged at 20,000 g and at 4°C for 30 min, and the supernatant was lyophilized. The lyophilisate was weighed to establish the recovery of LPS. LPS was reconstituted in water to make up an end
5 concentration of 5 mg/ml. Dependent of type of LPS and batch number, a volume of 1 mg/ml *porcine* hemoglobine (Hb) was added to a concentration as indicated in the text. Each batch was portioned into 0.5-mg LPS fractions, which were dried using a vacuum evaporator and then stored at 5-8°C.

Oxidation of LPS^[AAB7]

10 A portion of 0.5 mg hemoglobin-fortified LPS was dissolved in 500 µl 100 mM sodium acetate pH 5.5. Following the addition of 20 µl 50 mM sodium periodate (Sigma), the solution was incubated for 40 min on ice protected from light. The oxidation of LPS was quenched and the solution was desalted by passing 500 µl of the reaction mixture through a NAP-5 cartridge with a
15 gravity-controlled flow. Modified LPS was eluted with 1 ml 10 mM sodium acetate, pH 4.0. Prior to use, the cartridge was conditioned thrice with 3 ml 10 mM sodium acetate, pH 4.0.

Immobilization of LPS^[AAB8]

To immobilize the oxidized LPS antigens to the beads, the carboxylic groups at
20 the bead surface were activated with a mixture of EDC/NHS available from the amine-coupling kit for 20 min on a gyro rocker. Following centrifugation and removal of supernatant, activation was followed by a reaction with 5 mM aqueous carbohydrazide for 20 min. Beads with modified surface were pelleted again and upper liquid was discarded before addition of 1 M ethanolamine and
25 incubation for 20 min. Following another centrifugation step at 14,000 g for 5 min, oxidized LPS solved in sodium acetate pH 4.0 was added to allow immobilization for 90 min. Following removal of the suparenatant acquired through centrifugation, the linkage between bead-surface and antigen was

stabilized using 100 mM sodium cyanoborohydride solved in 10 mM sodium acetate at pH 4.

BioPlex assay

Following the warming-up of the bead counter device, this BioPlex (BioRad,
5 Veenendaal, The Netherlands) was calibrated according to the instructions of
the producer using a BioPlex calibration kit (BioRad). Samples were diluted in
50 mM PBS in wells of a microtiter plate which were then supplemented with
50 µL 5000 beads/mL LPS-coated beads. The antigen-antibody binding was
allowed for 30 min on a microtiter plate shaker operated at 200 rpm. Then
10 10 µL goat anti-swine IgG (H+L) tagged with Alexa532 fluorescent labels
diluted 8 times in 50 mM PBS were added and incubation was continued for 15
min on the shaker. Beads were then analysed for their fluorescence profiles
for 30 s on the BioPlex machine.

15 Results and discussion

Fluorescent beads were prepared for coating with LPS from different specific
Salmonella serovar sources representing serogroups B, C and D relevant as
zoonoses in foods from chicken and swine. It should be noted that serogroup E,
which is relevant for pork products, is not studied here. Following the
20 immobilization of each type of LPS to individual beads, which are internally
coded, the success of the coating was assessed using commercially available
monoclonal antisera against somatic antigens O4, O5, O7, O8 and O9.
However, while anti-O5 gave a response of 6398 units, anti-O9 gave 145 units,
whereas the background signal of non-matching antigens-antibodies was less
25 than 91 units in all cases (Fig. 6). In a similar way, anti-O4 and anti-O7 gave
responses of 305 units and 174 units, respectively (Fig. 7). These differences in
responses between commercially available antisera preparations were in very

good correspondence with those observed using a surface Plasmon resonance (SPR) biosensor and reflect differences in antibody titers.

In a similar way, the activity of identical LPS batches oxidized on different days were tested using a similar panel of commercial antisera (Fig. 8).

- 5 Compared to the other oxidation batch, responses ranged between 57% and 148%, which is susceptible for improvements.

Different preparations of meat drip, *i.e.* juice that is acquired from muscle tissue following a freeze and thaw cycle, and sera from chickens were analysed (Fig. 9). Recorded activities were as expected. Meat drip, serum and a
10 mixture of meat drip and serum from *Salmonella*-free chickens gave low abundant fluorescent conjugated beads. In contrast, anti-*S. pullorum* and anti-*S. gallinarum* should give a response on serogroups B and D, as it contains antigens O1 and O12, which it does for drip and serum. *S. infantis* contains antigen O6, which is shared by C₁ and C₂. Indeed, this activity is
15 observed in drip and serum[AAB9].

Besides chicken serum, prepared swine serum were tested as well (Fig. 10). Serum from *Salmonella*-free pigs gave MFI responses in the range of 110 units (serogroup C₂) to 137 units (serogroup B) and were close to the responses of beads only incubated with buffer, namely from 94 units (serogroup D) to 129
20 units (serogroup C₂). As expected, significant signals were recorded when sera were spiked with anti-*S. typhimurium* and *S. livingstone* antisera, namely 969 units on serogroup B and 207 units on serogroup C₁, respectively. The spiked sera did not react with non-corresponding antigens giving responses between 104 MFI units and 131 MFI units.

25

Example 5

Determination of anti-Salmonella antibodies in exotic avian species using an SPR biosensor[AAB10]

30

Aim of study

The aim of this investigation was to explore whether the developed SPR biosensor technology based on the use of immobilized selected *Salmonella* LPS to detect indirectly *Salmonella* infections in food-producing animals, is able to
5 detect such infections in exotic animal species as well.

Materials

Plasma from tocoucans (*Rhamphastos toco*) and a sharp-tailed grouse (*Tympanuchus phasianellus*), which were infected with *S. typhimurium* of
10 different phage types, were kindly provided by Dr. W. Schaftenaar and Ing. M. de Boer (Veterinary Department, Rotterdam Zoo, The Netherlands). The disease history of these animals is the following.

From the faeces of a Tocotoucan sampled at March 24, 1994, *S. typhimurium* phagetype 292 was isolated. From another faeces sample of the same bird,
15 *S. typhimurium* phagetype 352 was isolated at June 28, 1994. August 24, 1994, plasma was collected from this animal used for SPR analysis in this study.

Blood was collected from a diseased sharp-tailed grouse at October 28, 1997. The plasma prepared from this blood was used for analysis in this study. This
20 animal died the next day. *S. typhimurium* phagetype 507 was isolated from the dead bird.

Methods

An SPR biosensor (Biacore 3000) containing a sensor chip of which flow
25 channels were coated with LPS from *S. enteritidis*, *S. livingstone*, *S. goldcoast* and *S. typhimurium*, was operated as described earlier. Plasma samples were diluted as described for sera in examples 1 and 2 and analysed.

Results and discussion

The invention was in the first place developed for application in the food chains to secure the safety of food of animal origin with respect to *Salmonella* contaminations. Nonetheless, the field of applicability was tested with plasma collected from two exotic avian species, namely a tocotoucan (*R. toco*) and a sharp-tailed grouse (*T. phasianellus*), which were infected with *S. typhimurium* as disclosed by classic microbiological diagnostics (personal communication with M. de Boer, Blijdorp Zoo, Rotterdam).

The faeces of the tocotoucan was found positive five and two months before blood was sampled and a humoral response could develop over a relatively long period of time. Indeed, the biosensor response was high (Table 33). Compared to blank serum from SPF chickens and to standard antiserum (anti-serogroups A to S), the reactivity with *S. typhimurium* LPS was dramatic high (4185 response units (RU)). A response was also observed on the channel of *S. enteritidis* (1220 RU), which was also observed for serum from chickens highly infected with exclusively *S. typhimurium* and is in accordance with the presence of somatic antigen O12 in both serovars and thus in serogroups B and D (cf. Tables 1 and 3). Unexpectedly, a relatively high response was also observed at the *S. goldgoast* channel. It can not be excluded here that this bird was infected with multiple *Salmonella* serovars simultaneously or sequentially with *S. typhimurium* as the last infection, including a C₂ infection.

The plasma of the sharp-tailed grouse was not very reactive with the different LPS types, but reactivity was in all cases higher than that of the blank serum and on Sl and St LPS higher than that of the reference antiserum (Table 33). As the bird died rapidly from the infection, a significant humoral response against the selected antigens was probably not fully developed and not detected by the biosensor. It should be noted that serology is therefore not very suitable for diagnosis on an individual level. As evidenced by

Swanenburg (Utrecht Thesis, Utrecht University, 2000), serology is, in particular, suitable for assessing the Salmonella status on a population level.

Table 33. Results of the analysis of plasma collected from a tocotoucan and a sharp-tailed grouse which were infected with *S. typhimurium*. The results are expressed in relative response units. Samples were analysed thrice.

Flow channel coated with LPS from	Immobilis ation level of probing LPS	Serum from SPF ^a chickens	Anti- serogroup B ^b	Plasma from tocotouca n	Plasma from sharp- tailed grouse
Se	3241	17 (3)	-3 (1)	1220 (304)	33 (1)
Sg	3354	9 (4)	-8 (1)	492 (186)	59 (9)
Sl	5241	6 (6)	-29 (2)	58 (12)	89 (8)
St	2755	46 (6)	663 (33)	4185 (362)	73 (8)

^a serum from specific pathogen free chickens, i.e. blank serum

^b commercially available monoclonal antibody reacting with Salmonella serogroup B.

Example 6

Detection of bacteriophage Felix O1 (FO1) through its binding to Salmonella LPS immobilized on an SPR biosensor chip surface

5

Goal

To determine the binding of the bacteriophage FO1 to LPS immobilized to the biosensor surface

10 **Approach**

To prove binding of bacteriophage Felix O1 (FO1, Félix d'Hérelle Reference Centre for Bacterial Viruses, Laval, Canada) to *Salmonella* LPS, the bacteriophage was diluted in HBSEP to obtain a concentration series. These samples were injected for 2 min on the biosensor to allow binding to LPS from
15 *S. typhimurium*, *S. enteritidis*, *S. goldcoast* and *S. livingstone*, which were each immobilized separately on an individual flow channel of the sensorchip.

Results

The results are summarised in Fig. 11.

20

Although a relatively high concentration of bacteriophage was needed to obtain a significant response, it is evident from this experiment that 10^9 PFU of FO1 bacteriophages/mL and higher bound to LPS coupled to the chip surface.

25 **Example 7**

Detection bacteriophage FO1 through its binding to S. typhimurium LPS immobilized on an SPR biosensor following its incubation with S. typhimurium, S. enteritidis, S. goldcoast and S. livingstone

30

Goal

To determine the binding of bacteriophage FO1 to a live culture of *Salmonella* spp. in HBSEP.

Approach

- 5 Different concentrations of cultures of *S. typhimurium*, *S. enteritidis*, *S. goldcoast*, *S. livingstone* and blank medium were mixed for 5 min with 1.2×10^9 PFU bacteriophage FO1. After incubation, bacteria were spun down and supernatant was analysed on the Biacore with a biosensor chip containing immobilized LPS from *S. typhimurium*.

10

Results

- In order to determine significant responses, a cut off value was established by the averaged readings of blank medium containing no *Salmonella* but 1.2×10^9 PFU bacteriophage, minus 3 times standard deviation. Applying this value
15 disclosed that *Salmonella* should be present at a concentration of at least 6×10^8 CFU/mL, 3×10^6 CFU/mL, 4×10^7 CFU/mL and 3×10^4 CFU/ml for *S. typhimurium*, *S. enteritidis*, *S. goldcoast*, *S. livingstone*, respectively, to give a significant response (Figure 12).

20 Discussion

- The absorption rate of bacteriophage FO1 to *Salmonella* spp. is probably dependent of the accessibility of *N*-acetylglucosamine in the core region, the binding site of FO1 (Lindberg, 1977; Lindberg and Holme, 1969[AAB11]). Long and numerous O-side chains occurring in the polysaccharide region of the LPS
25 of targeted *Salmonellae* could, therefore, impair the binding of FO1 to the analyte. It is, therefore, expected that free bacteriophages will have a variable binding capacity towards *Salmonella* strains and that propagation of the virus will largely depend on the molecular profile of the exposed LPS.

- When immobilizing LPS-protein complexes at the surface of a solid carrier for
30 a diagnostic method, as described in the invention, a dense network of ligands

may be formed. In the presented Biacore analysis, the density of the complex and the hindrance by the proteins may play a role in the observed difference in bacteriophage binding to the four LPS types.

It should be noted here that, as discussed in the invention as well, oxidation of
5 monosaccharide constituents in the core region is expected, including that of the N-acetylglucosamine. The ligand for the bacteriophage is, therefore, affected and this may influence the sensitivity of the test.

10 Example_[AAB12] 8

Propagation of bacteriophage FO1 in Salmonella and non-Salmonella strains

Goal

15 The propagation of FO1 is not exclusive in *Salmonella* spp., but possibly also in non-*Salmonella* strains (3). This study was initiated to investigate the selectivity of the of proliferation of the bacteriophage FO1. For this purpose, a number of important food non-*Salmonella* pathogens were exposed to the FO1 bacteriophage.

20

Methods

Seven non-*Salmonella* strains (*Listeria monocytogenes*, *Escherichia coli*,
Pseudomonas aeruginosa, *Bacillus cereus*, *Citrobacter*, *Enterococcus faecalis*,
and *Staphylococcus aureus*) and seven *Salmonella* strains (*S. cholerasuis*,
25 *S. berta*, *S. meleagridis*, *S. Agona*, *S. pullorum*, *S. Virchow* and *S. enteritidis*)
were grown in Trypton Soy Broth (Oxoid CM129).

At t 0 hours, 1.2×10^8 PFU of FO1 (end concentration 1×10^6 PFU/ml) was added to all cultures. Every hour a sample was drawn, and absorbance at λ 600 nm (Figs. 13 and 14), plaque forming units/ml (Figs.15 and 16) and, after

concentration and buffer exchange, binding to immobilized St-LPS immobilized on a sensorchip surface in a Biacore biosensor (Fig. 17), were determined.

- 5 Some non-*Salmonella* bacteria, including *L. monocytogenes*, *P. aeruginosa*, *E. faecalis* and *Staph. aureus*, did not show growth in five hours of culture (Fig. 13). These bacteria were obviously not loaded nor infected by the bacteriophage FO1, because the concentration of bacteriophages did not rise over time and was stationary at 1×10^6 PFU/ml (Fig. 15).
- 10 In contrast, all *Salmonella* serovars, except *S. virchow*, grew in five hours of incubation (Fig. 14). This *S. virchow* strain was probably lysed completely by the proliferating bacteriophages, as a clear increase in the concentration of bacteriophages can be shown from 1×10^6 PFU/ml to 1×10^{10} PFU/ml (Fig. 16). In a similar way, *S. berta* and *S. meleagridis* bacteria showed an increase of
- 15 bacteriophages concentration. However, these bacteria, in particular, *S. berta* showed good growth (Fig. 14).

In the context of the invention, the binding of propagated bacteriophages to the sensor surface is of greatest interest. For this purpose, bacteriophages propagated in *Salmonella* were concentrated, dialysed and serially diluted

20 before SPR biosensor analysis (Fig. 17). Unexpectedly, comparable bacteriophage concentrations gave dramatic different biosensor responses. Most probably phages were lost, in particular during concentration and dialysis step, during the sample preparation. Nevertheless, the preparations of bacteriophages, which have shown propagation (cf. Fig. 16) gave clearly

25 higher responses than the blank sample, which contained the starting concentration of bacteriophages only.

Conclusion

These experiments showed that bacteriophages can be used as an analytical tool to detect the presence of *Salmonella* in samples and that *Salmonella* LPS immobilized to a solid surface can be used to probe the increment of the bacteriophages as a result of propagation of the virus in the host bacteria
5 following a short incubation period.

Example 9

*Immobilisation of salmonella-derived LPS onto fluorescent beads (luminex)
10 and detection of antibodies reporting a current or past salmonella infection in various biological samples*[AAB13]

0. WARNING AND SAFETY PRECAUTIONS

Lipid polysaccharides are potent immunogens, which can bring sensitive
15 persons into a septic shock upon intake or inhalation. Precautions should be made to prevent contact with this biochemical.

Sodium periodate is an oxidizing agent and may cause explosions when brought in contact with strong reducing agents.

The procedure utilizes sodium cyanoborohydride. The procedure should
20 therefore be carried out with precautions as with for example hand gloves and a mask. Use this chemical substance only in a chemical fume hood. The material is very toxic to aquatic organisms and may cause long-term adverse effect in the aquatic environment. The material and solution waste should be disposed of as hazardous waste.

25

1. INTRODUCTION

Among the components of the antigenic structure of the genus *Salmonella*, the somatic antigens are important as an instrument to trace immune response in

animals upon an invasive infection of this micro-organism. Somatic antigens are located on the polysaccharide part of lipid polysaccharide (LPS), which is a constituent of the bacterial cell wall. After extraction and isolation of LPS from carefully chosen *Salmonella* serotypes (*cf.* SOP CHEMIE/A21[AAB14]), antigen-
5 containing LPS is coupled covalently to beads, which are internally coded by a specific mixture of fluorescent material. The exploited technology platform from Luminex can identify up to 100 differently internally coded beads in a single test. A single species of beads can be immobilized with a mixture of LPS reflecting different serogroups, or different species of beads can be each
10 immobilized with LPS reflecting a single specific serogroup or serovar of the pathogenic micro-organism. The LPS-containing beads are incubated with body-derived materials, such as blood, plasma, serum, meat drip/juice, egg-yolk, milk etc , to enable anti-*Salmonella* antibody-antigen binding. The specific binding is detected following a second incubation with fluorescently
15 tagged anti-immunoglobulin antibodies in a device analyzing simultaneously the emission wavelengths of excited beads and tagged antibodies. Only when both fluorescence of bead and antibody are detected simultaneously a response to a specific bead species will be recorded.

This SOP describes the method for oxidation, immobilization of LPS onto
20 beads (Luminex) and an assay to assess the quality of produced.

2. SCOPE AND FIELD OF APPLICATION

To analyze biological fluids, such as serum samples from chicken and pigs, for the presence of anti-*Salmonella* antibodies reacting with O3, O4, O5, O6, O7,
25 O8, O9, O10 and O12 somatic antigens reflecting a history of or current infection of *Salmonella* from serogroups B, C, D and E.

3. REFERENCES

Extraction and isolation of lipopolysaccharides from *Salmonella* spp.;

Immobilisation of salmonella-derived LPS onto a biosensor chip (biacore) and

5 detection of serum antibodies reporting a current or past salmonella infection;

Optimalisation of protein addition to LPS for immobilization and detection of serum antibodies; See example 1.

4. DEFINITIONS

10 c = concentration in % (m/v), % (v/v), mol/l or mmol/l as indicated.

5. LIST OF ABBREVIATIONS

5.1 EDC, 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide hydrochloride

5.2 LPS, Lipopolysaccharides

15 5.3 NaCl, Sodium chloride

5.4 NaOH, Sodium hydroxide

5.5 NHS, *N*-hydroxysuccinimide

5.6 PBS, Phosphate buffered saline

5.7 Se, *Salmonella enteritidis*

20 5.8 Sg, *Salmonella goldcoast*

5.9 Sl, *Salmonella livingstone*

5.10 Sm, *Salmonella meleagridis*

5.11 SOP, Standard operating procedure

5.12 St, *Salmonella typhimurium*

25 5.13 TCA, Trichloric acetic acid

6. PRINCIPLE

LPS is oxidized in the presence of a protein facilitated by sodium periodate. The LPS-protein solution is desalted using a NAP-5 column. After activation of
5 the carboxylic acid groups at the surface of the beads with the aid of 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide hydrochloride (EDC),
N-hydroxysuccinimide (NHS) followed by a reaction with carbonylhydrazide, desalted oxidized LPS-protein complex is immobilized to the solid phase of the beads. Bound LPS is then stabilized with sodium cyanoborohydride. Prior to
10 routinely use, the performance of bead-conjugated LPS to bind anti-*Salmonella* antibodies is assessed using a panel of reference monoclonal agglutination sera.

7. REACTIONS

15 7.1 OXIDATION OF CARBOHYDRATE MOIETY

Periodate will induce an oxidative disruption of linkages between vicinal cis-diols on, in particular, carbohydrate moieties, as in *e.g.* mannose, to yield aldehyde functionalities, see Figure 18. This reaction is typically performed in buffers at a pH range between 4.5 and 5.5 in the dark using a freshly prepared
20 10-100 mM sodium meta-periodate in 0.1 M sodium acetate.

NOTE 1: The positions of conjugations indicated in Figure 18, to link up the depicted monosaccharide with other monosaccharide residues in a polysaccharide, as in *e.g.* LPS, are here just given as an example. **R'** and **R**
25 indicate the distal and the proximal positions, respectively, in the carbohydrate chain. The oxidation of hydroxyl groups into aldehydes may repeat itself within the polysaccharide chain in each monosaccharide constituent containing susceptible vicinal diols.

NOTE 2: Periodate will also oxidize, when present, certain β -aminoethanol derivatives such as the hydroxylysine residues in collagen, as well as methionine (to its sulfoxide) and certain thiols (usually to disulfides). In addition, N-terminal serine and threonine residues of peptides and proteins
5 can be selectively oxidized by periodate to aldehyde groups. These reactions, however, usually occur at a slower rate than oxidation of vicinal diols.

7.2 CONJUGATION TO PROTEIN

Oxidation is performed in the presence of a protein. The bis-aldehyde
10 compounds, such as the oxidized monosaccharide constituents in the polysaccharide chain of LPS here, may react with any amino group in a protein and may form a Schiff-base linkage resulting in a substituted imine. See Figure 19.

NOTE: The substituted imine is stabilized while the complex is attached to
15 the bead surface, by a reduction facilitated by cyanoborohydride. This type of reaction scheme is known as a reductive amination.

7.3 IMMOBILIZATION TO FLOURESCENT BEADS

The carboxylic acid labeled beads are activated using *N*-ethyl-*N*-(3-dimethyl
20 aminopropyl)-carbodiimide hydrochloride (EDC) and *N*-hydroxysuccinimide (NHS). The activation is followed by a reaction with carbohydrazine. The reactive aldehyde functionalities react spontaneously with the hydrazide to hydrazones, which are then reduced to stabilise the covalent bonds. See Figure
20.

25

NOTE: The protein (R'') carries multiple $-NH_2$ groups and can therefore be conjugated with multiple oxidized LPS entities. At the other hand, the polysaccharide part in LPS may carry multiple free aldehyde groups in a

single molecule. These aldehyde groups may for a part or completely captured by the hydrazide-layer on the beads. The net result may be a very stable complex network of protein-LPS covalently linked to the bead surface.

5 8. REAGENTS AND MATERIALS

In the complete procedure only reagents of recognized analytical grade and only distilled water or water of equivalent purity are used, unless stated otherwise. Reference to a company is for information and identification only and does not imply a recommendation unless so stated.

10 8.1 CHEMICALS

8.1.1 Acetic acid (J.T. Baker, Deventer, The Netherlands)

8.1.2 Amine coupling kit (Biacore AB, Uppsala, Sweden) consisting of:

8.1.2.1. Vial containing 115 mg *N*-hydroxysuccinimide (NHS)

8.1.2.2. Vial containing 750 mg 1-ethyl-3-(3-dimethylaminopropyl)

15 carbodiimide hydrochloride (EDC)

8.1.2.3. Vial containing 10.5 ml, $c = 1$ mol/l, ethanolamine hydrochloride – sodium hydroxide pH 8.5

8.1.3. Bio-Plex Calibration kit (Bio-Rad, Veenendaal, the Netherlands)

8.1.4. Carbohydrazide, $\text{CN}_4\text{H}_6\text{O}$ (Fluka Chemie GmbH, Buchs, Switzerland)

20 8.1.5. Carboxymethyl-dextran sodium salt (Fluka)

8.1.6. Potassium dihydrogen phosphate (KH_2PO_4) (Merck, Darmstadt, Germany)

8.1.7. Proclin 150 (Supleco, Bellefonte, PA, USA)

8.1.8. Monoclonal anti *Salmonella* O-antigens:

25 8.1.8.1. anti-O4 (SIFIN, Berlin Germany)

8.1.8.2. anti-O5 (SIFIN)

8.1.8.3. anti-O6₁ (SIFIN)

8.1.8.4. anti-O7 (SIFIN)

8.1.8.5. anti-O8 (SIFIN)

8.1.8.6. anti-O9 (SIFIN)

5 8.1.8.7. anti-O10 (SIFIN)

8.1.9. *Salmonella* LPS, in-house isolated LPS by TCA extraction (SOP CHEMIE/A21) prepared from the *Salmonella* bacteria serovars *enteritidis* (Se), *goldcoast* (Sg), *livingstone* (Sl), *meleagridis* (Sm) and *typhimurium* (St) with protein (SOP CHEMIE/A23)

10 8.1.10. Sheep anti-mouse Ig-PE (Chemicon, Boronia, Victoria, Australia)

8.1.11. Sodium acetate trihydrate (J.T. Baker, Phillipsburgh, NJ, USA)

8.1.12. Sodium chloride (Merck)

8.1.13. Sodium cyanoborohydride (NaCNBH₃) (Fluka)

8.1.14. di-Sodium hydrogen phosphate (Na₂HPO₄) (Merck)

15 8.1.15. Sodium hydroxide, *c* = 50 mmol/l (Biacore)

8.1.16. Sodium *m*-periodate (NaIO₄) (Sigma-Aldrich, Zwijndrecht, the Netherlands)

8.1.17. Water is obtained from a Milli Q water purification system

20 8.2. SOLUTIONS

8.2.1. Acetic acid solution, *c* = 0.1 g/ml: Dilute 1 ml acetic acid (8.1.1) with 9 ml water.

8.2.2. Acetate buffer solution, *c* = 10 mmol/l, pH 4.0: dissolve 0.272 g sodium acetate trihydrate (8.1.11) in 180 ml water and adjust to pH 4.0 with acetic acid (8.1.1) and make up to 200 ml with water. This buffer is stable for approximately six months.

25

- 8.2.3. Acetate buffer solution, $c = 1.0 \text{ mol/l}$, pH 5.5: dissolve 13.6 g sodium acetate trihydrate (8.1.11) in 90 ml water and adjust to pH 5.5 with acetic acid (8.1.1) and make up to 100 ml with water. This buffer is stable for approximately six months.
- 5 8.2.4. Acetate buffer solution, $c = 100 \text{ mmol/l}$, pH 5.5: dilute 1.0 ml acetate buffer solution, $c = 1.0 \text{ mol/l}$ (8.2.3) with 9.0 ml water.
- 8.2.5. Carbohydrazide solution, $c = 100 \text{ mmol/l}$: dissolve 9.0 mg carbohydrazide (8.1.4) in 1000 μl water.
- 8.2.6. Carbohydrazide solution, $c = 5 \text{ mmol/l}$: dilute 10 μl carbohydrazide solution (8.2.5) with 190 μl water. Prepare just before use.
- 10 8.2.7. EDC-solution: reconstitute EDC (8.1.2.2.) in 10.0 ml water.
- 8.2.7.1. Fractionate 100- μl aliquots of this solution (8.2.7) in polypropylene tube (9.15). Store at -18°C or at lower temperature. The aliquots are stable for two months. Before use: Thaw frozen aliquots and agitate them gently to ensure homogeneous solutions.
- 15 8.2.8. Ethanolamine solution: Pipette 200 μl $c = 1 \text{ mol/l}$ ethanolamine solution (8.1.2.3) in a polypropylene tube (9.15)
- 8.2.9. NHS-solution: reconstitute NHS (8.1.2.1) in 10.0 ml water.
- 8.2.9.1. Fractionate 100- μl aliquots of this solution (8.2.9) in polypropylene tube (9.15). Store at -18°C or at lower temperature. The aliquots are stable for two months. Before use: Thaw frozen aliquots and agitate them gently to ensure homogeneous solutions.
- 20 8.2.10. PBS (5.6), $c = 100 \text{ mmol/l}$, pH 7.2: dissolve 67.9 g sodium chloride (8.1.12), 14.7 g di-sodium hydrogen phosphate (8.1.14) and 4.3 g potassium dihydrogen phosphate (8.1.6) in 1 L water (8.1.17).
- 25 8.2.11. PBS, $c = 10 \text{ mmol/l}$, pH 7.2: dilute 100 mL 10 times concentrated PBS (8.2.10) in 1 L water (8.1.17).

- 8.2.12. Anti-mouse Ig-PE, prediluted: dilute fluorescent conjugate 5 times by mixing 40 μ l anti-mouse Ig-PE (0) with 160 μ l PBS (8.2.11).
- 8.2.13. Sodium cyanoborohydride, $c = 1.00$ mol/l: dissolve 62.8 mg sodium cyanoborohydride (8.1.13) in 1000 μ l 10 mM acetate solution pH 4.0 (8.2.2).
- 8.2.14. Sodium cyanoborohydride, $c = 100$ mmol/l: dilute 180 μ l sodium cyanoborohydride solution (8.2.13) with 1620 μ l acetate solution, $c = 10$ mmol/l, pH 4.0 (8.2.2). Prepare just before use.
- 8.2.15. Sodium hydroxide, $c = 5$ mmol/l: Dilute 400 μ l sodium hydroxide $c = 50$ mmol/l (8.1.15) with 3.6 ml water in glass vial (9.10).
- 8.2.16. Sodium *m*-periodate solution, $c = 100$ mmol/l: Add to 214 mg sodium *m*-periodate (8.1.16) 10.0 ml water (8.1.17).
- 8.2.17. Sodium *m*-periodate 'ready to use': Pipette 100 μ l of sodium *m*-periodate solution, $c = 100$ mmol/l (8.2.16) in a 1.4 ml polypropylene tube (9.14) and dry with a centrifugal evaporator (9.4).
- 8.2.18. Sodium periodate solution, $c = 50$ mmol/l: Dissolve sodium periodate 'ready to use' (8.2.17) in 200 μ l acetate solution, $c = 100$ mmol/l pH 5.5 (8.2.4). Prepare just before use.

20 8.3. STANDARD MONOCLONAL REFERENCE SOLUTION

- 8.3.1. Concentrated monoclonal *Salmonella* anti-O5 (8.1.8.2): dilute 5 μ l anti-O5 10 times by adding 45 μ l PBS $c = 10$ mmol/l, pH 7.2 (8.2.11).
- 8.3.2. Monoclonal *Salmonella* anti-O5: 7.5 μ l anti-O5 (8.3.1) 10 times diluted by addition of 67.5 μ l PBS $c = 10$ mmol/l, pH 7.2 (8.2.11).
- 25 8.3.3. Monoclonal anti *Salmonella* O-antigens (8.1.8): dilute 7.5 μ l of each monoclonal (8.1.8.1, 8.1.8.3, 8.1.8.4, 8.1.8.5, 8.1.8.6, 8.1.8.7) with 67.5 μ l PBS (8.2.11) in a micro titerplate (9.20).

8.3.4. Thrice diluted monoclonal antibodies: add 25 µl monoclonal antibody solution (8.3.2 and 8.3.3) to 50 µl PBS $c = 10$ mmol/l, pH 7.2 (8.2.11) in wells of the same micro titerplate (8.3.3)

5 8.3.5. Repeat step 8.3.4 twice to obtain 9 and 27 times diluted antibodies in fresh wells of the micro titerplate in the case of anti-O4, anti-O6, anti-O7, anti-O8, anti-O9 and anti-O10. In the case of anti-O5, these dilution factors were 90 and 270 times, respectively.

8.3.6 Remove 25 µl from the highest dilution (8.3.5).

8.3.7 The antibody solutions are now ready for use.

10

NOTE The final dilution factors are 100, 300, 900 and 2700 times in the case of anti-O5, whereas in the other cases antibodies were finally diluted 10, 30, 90 and 270 times compared to the original preparation (8.1.8).

15 8.4. AUXILIARY MATERIALS

8.4.1 NAP-5 column (0.5 ml, Sephadex G-25, Amersham Biosciences).

8.4.2. COOH-Beads (5.6 µm COOH microspheres) numbers 24, 25, 26, 27 and 28 mixed in 0.01% aqueous merthiolate at 1.25×10^7 beads/mL (BioRad).

20 9. EQUIPMENT AND APPARATUS

Reference to a company is for information and identification only and does not imply a recommendation unless so stated. Equivalent equipment and apparatuses may be as appropriate as well.

9.1. Analytical balance (type AE 240, Mettler, Zürich, Switzerland)

25 9.2. BioPlex 200 (Bio-Rad)

9.3. Bürker-Türk counting chamber (Mainit B.V., Wormer, The Netherlands)

- 9.4. Centrifugal evaporator (Jouan, Saint-Herblain, France)
- 9.5. Finn pipette, 5-40 μ l (Labsystems Oy, Helsinki, Finland)
- 9.6. Finn pipette, 40-200 μ l (Labsystems Oy)
- 9.7. Finn pipette, 200-1000 μ l (Labsystems Oy)
- 5 9.8. Finn pipette, 1- 5 ml (Labsystems Oy)
- 9.9. Glass collection tube, 5 ml, stoppered (Renes, Zeist, The Netherlands)
- 9.10. Glass vial diameter 16 mm (Sarstedt B.V., Etten-Leur, the Netherlands)
- 9.11. Gyro rocker SSL3 (Stuart, Barloworld Scientific Ltd. Stone,
10 Staffordshire, United Kingdom)
- 9.12. Microscope Axiovert 25 (Zeiss, Sliedrecht, the Netherlands)
- 9.13. Microscope glass cover slides (Chance Propper Ltd, Smethwick, Warley, England)
- 9.14. Mini shaker with microtiter cup head (MS 1, Janke & Kunkel, Staufen,
15 Germany)
- 9.15. Polypropylene tube, diameter 7 mm, with snap caps (Biacore)
- 9.16. Polypropylene tube, 1.4 ml (Micronic, Lelystad, The Netherlands)
- 9.17. pH meter (type pH 537, WTW, Weilheim, Germany).
- 9.18. Sonification bath, Branson 2200 (Branson Ultrasonics B.V., Soest, The
20 Netherlands)
- 9.19. Vacuum manifold, to run several NAP-5 columns simultaneously (type SPE-12G, with PTFE stopcock(s), J.T.Baker).
- 9.20. V-bottomed microtiter polystyrene plate, 96-well format (Greiner Bio-one GmbH, Frickenhausen, Germany)
- 25 9.21. Vortex-mixer (KS-1, Janke & Kunkel)

10. SOFTWARE

The BioPlex apparatus is operated with Bio-Plex Manager software 4.1.

11. PROCEDURE

5 11.1. OXIDATION AND DESALTING OF LPS SOLUTION

SAFETY PRECAUTION:

Lipid polysaccharides are potent immunogens, which can bring sensitive persons into a septic shock upon intake or inhalation. Precautions should be made to prevent contact with this biochemical.

10 Sodium periodate is an oxidizing agent and may cause explosions when brought in contact with strong reducing agents.

The procedure utilizes sodium cyanoborohydride. The procedure should therefore be carried out with precautions, such as using hand gloves and a mask. Use this substance only in a chemical fume hood.

15 The material is very toxic to aquatic organisms and may cause long-term adverse effect in the aquatic environment. The material and solution waste should be disposed of as hazardous waste.

11.1.1. OXIDATION

11.1.1.1. Add 500 µl acetate buffer $c = 100$ mmol/l, pH 5.5 (8.2.4) to dry LPS
20 (8.1.9; See safety precaution)

11.1.1.2. Vortex (9.21) the solution (11.1.1.1) and sonicate (9.18) for 20 min and observe the reconstitution process so that all LPS is dissolved.

11.1.1.3. Add 20 µl periodate solution $c = 50$ mmol/l (8.2.18) to the LPS solution (11.1.1.2)

25 11.1.1.4. Vortex (9.21) the solution (11.1.1.3)

11.1.1.5. Incubate on ice for 40 min protected from light.

11.1.1.6 Quench oxidation by desalting the solution (11.1.1.4) as described in

11.1.2

11.1.2. DESALTING

5 11.1.2.1. Place NAP-5 column(s) (8.4.1) on manifold (9.19).

11.1.2.2. Condition the column(s) (11.1.2.1) by passing three 3-ml portions of acetate buffer $c = 10$ mmol/l, pH 4.0 (8.2.2) over the column bed on a flow generated by gravity only. Allow the buffer to enter the gel bed completely.

10 11.1.2.3. Pipette 0.5 ml oxidized LPS solution (11.1.1.5) on the column. Allow the sample to enter the gel bed completely. The flow-through is not collected.

11.1.2.4. Elute oxidized LPS with 1 ml acetate buffer $c = 10$ mmol/l, pH 4.0 (8.2.2). Collect eluate in a 5-ml glass tube (9.9).

15 11.1.2.5. Vortex (9.21) the solution (11.1.2.4) for 10 s and add 2 μ L Proclin 150 (8.1.7).

11.1.2.6. When not immediately used (11.1.2.5) store samples at 4°C to 7°C.

11.1.2.7. Prior to immobilization, the LPS-containing solution (11.1.2.5) which can be used for different matrix and species applications is diluted as
20 indicated in the following Tables 34 and 35.

Table 34. Amount of LPS solution used to immobilize beads for detection of antibodies to Salmonella O-antigens in swine sera

LPS type (8.1.9)	LPS stock solution (11.1.2.6) in μl	Added volume of acetate buffer, pH 4.0 (8.2.2) in μl
Se	75	225
Sg	75	225
Sl	75	225
Sm	75	225
St	75	225

Table 35. Amount of LPS solution used for the immobilization to fluorescent beads for detection of antibodies in chicken sera reacting with Salmonella O-antigens.

LPS type (8.1.9)	LPS stock solution (11.1.2.6) in μl	Added volume of acetate buffer, pH 4.0 (8.2.2) in μl
Se	150	150
Sg	150	150
Sl	150	150
Sm	150	150
St	150	150

11.2. IMMOBILIZATION OF LPS TO BEADS

11.2.1. Beads (8.4.2) are vortex-mixed for minimally 1 min

10 11.2.2. Transfer a portion of 300 μ L beads (11.2.1) into a fresh container (9.16)

- 11.2.3. Centrifuge at 14,000 *g* for 5 min
- 11.2.4. Remove supernatant carefully from the beads using a 200 μ l pipette (9.6).
- 11.2.5. Leave a small amount of solution (10 μ L) in the vial and mix the beads
5 in the remaining solution on a vortex (9.21).
- 11.2.6. Thaw two portions of 100 μ l EDC (8.2.7.1)
- 11.2.7. Thaw two portions of 100 μ l NHS solution (8.2.9.1).
- 11.2.8 Mix 180 μ l of EDC (11.2.6) and 180 μ l of NHS (11.2.7).
- 11.2.9. Transfer 300 μ L EDC/NHS mix (11.2.8) to the beads (11.2.4) and
10 suspend rigorously using a pipette.
- 11.2.10. Facilitate reaction on a gyro rocker (9.11) for 20 min.
- 11.2.11. Centrifuge at 14,000 *g* for 5 min.
- 11.2.12. Carefully remove supernatant from beads, leave a small volume
(ca. 10 μ L) on top of pellet and suspend beads using vortex mixer (9.21).
- 15 11.2.13. Add 300 μ L 5 mM carbonyldiimidazole solution (8.2.6) to the beads
(11.2.12) and suspend rigorously using a pipette.
- 11.2.14. Facilitate reaction on a gyro rocker (9.11) for 20 min.
- 11.2.15. Centrifuge at 14,000 *g* for 5 min, remove supernatant from beads,
leave a small volume (ca. 10 μ L) on top of pellet and suspend beads
20 using vortex mixer (9.21).
- 11.2.16. Add 300 μ L 1 M ethanolamine solution (8.2.8) to the beads (11.2.15)
and suspend rigorously using a pipette.
- 11.2.17. Facilitate reaction on a gyro rocker (9.11) for 20 min.
- 11.2.18. Centrifuge at 14,000 *g* for 5 min, remove supernatant from beads,
25 leave a small volume (ca. 10 μ L) on top of pellet and suspend beads

- using vortex mixer (9.21).
- 11.2.19. Add 300 μ L diluted oxidized LPS in sodium acetate $c = 10$ mmol/l, pH 4.0 (11.1.2.7) and suspend rigorously using a pipette.
- 11.2.20. Allow reaction on a gyro rocker (9.11) for 90 min.
- 5 11.2.21. Centrifuge at 14,000 g for 5 min, remove supernatant from beads, leave a small volume (ca. 10 μ L) on top of pellet and suspend beads using vortex mixer (9.21).
- 11.2.22. Add 300 μ L cyanoborohydride solution $c = 100$ mmol/l (8.2.14) and suspend rigorously using a pipette.
- 10 11.2.23 Facilitate reaction on a gyro rocker (9.11) for 60 min.
- 11.2.24. Centrifuge at 14,000 g for 5 min and carefully remove supernatant from beads, leave a small volume (ca. 10 μ L) on top of pellet and suspend beads using vortex mixer (9.19).
- 11.2.25. Add 300 μ L PBS (8.2.11).
- 15 11.2.26. Add 1 μ L Proclin 150 (8.1.7) and mix suspension
- 11.2.27 The beads are ready for testing *Salmonella* antibodies in biological materials
- 11.2.28. Counting of LPS coupled beads
- 11.2.28.1. Vortex LPS coupled beads suspensions (11.2.25) and transfer
- 20 1 μ L in a fresh 1-mL tube (9.15)
- 11.2.28.2. Dilute by adding 24 μ L PBS $c = 10$ mmol/l, pH 7.2 (8.2.11) and mix.
- 11.2.28.3 The external supports of a Bürker-Türk (9.3) counting chamber are to be moistened with milliQ water (8.1.17) and the cover glass (9.13) is gently pushed onto the counting chamber from the front.
- 25 11.2.28.4 Fill a pipette with 20 μ L bead solution (11.2.28.2), gently form a drop at the tip of the pipette.
- 11.2.28.5 This drop (11.2.28.4) is to be placed between the cover glass and the

counting chamber.

- 5 11.2.28.6 As a result of the capillary effect the gap between the cover glass and the chamber base fills up. Before the bead solution can overflow at the edges of the chamber section, the tip of the pipette must be removed. If any air bubbles are visible or if the liquid has overflowed over the edges and into the grooves, the chamber must be cleaned and feeding must be repeated.
- 11.2.28.7 Place the filled counting chamber under a microscope (9.12) and magnify the image with a 10x object.
- 10 11.2.28.8 The count should be started at the top left-hand corner and follow the direction shown by the arrow (Figure 23, lower panel). Counting may be enhanced with the microscopes illumination reduced.
- 11.2.28.9 Count the number of beads in 16 squares (Figure 23, upper panel) inside the thick lined area (see Figure 24).
- 15 11.2.28.10 Notes on counting:
- 11.2.28.10.1 Use reduced microscope illumination for all chambers.
- 11.2.28.10.2 The difference of the counter cells in the large squares and the group squares must not exceed 10 cells.
- 20 11.2.28.10.3 Double checks must be performed for all cell counts. After counting the two counting nets the bottom counting net is to be counted in the same way as a check. When doing this it is to be ensured that the chamber has not dried out. This can be prevented by filling the bottom chamber only shortly before the count and the counting after the sedimentation time.
- 25 11.2.28.10.4 The difference between the totals of the counts for the two counting nets must not exceed 10 cells. The average value of the counts is then used in the calculation formula or multiplied by the corresponding factor.
- 30 11.2.28.11 Multiply the counted number (11.2.28.9) with the dilution factor (25 x) divided by counted area (1 mm²) multiplied with chamber depth to calculate the concentration of beads per ml.

11.3. DETECTION OF ANTI-SALMONELLA ANTIBODIES

Note: Ensure that all solvent and reagent reservoirs contain sufficient volume for complete method run before initiating the assay.

- 5 11.3.1. Make the BioPlex apparatus (9.2) operational by a 30 minutes laser warming up step, followed by a start up and calibration procedure with appropriate calibration solvents (8.1.3) according to the quick guide.
- 10 11.3.2. Mix and dilute LPS coated beads (11.2.26) with PBS $c = 10 \text{ mmol/l}$, pH 7.2 (8.2.11) so that the concentration of each bead equals 5000 per ml.
- 11.3.3. Transfer 50 μL bead mix (11.3.2) into a micotiter plate already filled with diluted monoclonal anti-O antigens (8.3.7).
- 11.3.4. Incubate for 30 min on a microtiter plate shaker (9.14).
- 15 11.3.5. Add 10 μL 5 times diluted anti-mouse Ig-PE (8.2.12).
- 11.3.6. Incubate for 15 min on a microtiter plate shaker (9.14).
- 11.3.7. Note in logging sample wells and their contents.
- 11.3.8. Place plate in BioPlex (9.2), set maximal counting time to 120 s and count at least 50 beads per LPS group.
- 20 11.3.9. Activate software program to count fluorescence of beads, which had captured (fluorescent) antibodies.

For a schematic representation of the procedure see Figure 21.

Typical responses of salmonella monoclonal antibodies are presented in Table 36.

Table 36. Typical responses of reference sera incubated with LPS-coated beads and a secondary fluorescent antibody providing the signal.

	LPS B	LPS C ₁	LPS C ₂	LPS D
anti-O4	305	76	95	88
anti-O5	6398	78	91	79
anti-O7	86	174	91	98
anti-O8	90	86	1668	91
anti-O9	82	81	82	145

Example 10

MILD PERIODATE OXIDATION

The success of the final binding of anti-*Salmonella* antibodies, und thus the
5 screening of invasive infections in the animal, is much dependent of the
oxidation of the monosaccharide constituents of the polysaccharide part of
LPS, and the oligossacharide part of the core region of LPS. It can be deduced
from the described structures for the different serotypes of e.g. *Salmonella* that
oxidation may lead to a breakdown of the antigenic structures, which are, in
10 particular, part of the polysaccharide part of LPS.

Whereas the oxidation of hexitols occurs rapidly, parynosides, which are
predominantly occurring in *Salmonella* LPS, need a higher periodate
concentration to facilitate the oxidation in the same time. Pyranosides, which
possess α -erythro-hydroxyl groups, such as in *arabino*, *galacto* or *manno*
15 configurations like in *Salmonella* spp. LPS, are easier oxidized than α -threo-
hydroxyl groups, such as in *xylo* or *gluco* variants. It should be realized that
while the ring is opened and aldehyde functions for attachment of e.g. protein
molecules are created, also α -hydroxy carbonyl compounds may be created,
which may oxidize again if periodate is still present. It is, therefore, that non-
20 conjugated monosaccharides, which are thus not part of an oligo- or
polysaccharide, are completely destroyed by a periodate oxidation to formic
acid and formaldehyde at sufficient high concentrations of the oxidizer. At
relatively high concentrations of periodate, 1,3-diketones and also di-axial
diols can be oxidized.

25 A breakdown or oxidation reaction more than only the creation of aldehyde
groups would, therefore, lead to a failure to detect an infection in a sample of
biological material despite a good coupling reaction to a solid phase supported
by the presence of a polyamine-containing molecule, such as a protein. In

other words, a mild periodate reaction is needed to leave the antigenic structure intact, but just enough to enable a coupling between protein and thus solid phase.

5 RESULTS:

LPS from Se, Sg, Sl and St was oxidized for 40 min at pH 5.5 using a range of sodium *m*-periodate concentrations. Following oxidation, LPS was coupled to a biosensor surface and immobilization efficiency and antigenic activities was monitored.

10

From Fig. 25 it is obvious that a higher oxidation grade of LPS from *S. enteritidis* gave rise to a corresponding higher coating level at the biosensor chip. However, despite the higher immobilization levels, the response from the antibody probing decreases as demonstrated in Figs. 26 and 27. For this LPS
15 type an optimum periodate concentration of 1.8 mM was determined.

In a similar way, the effects of oxidation on the immobilization and antigenic activity of LPS from *S. goldcoast* were tested (Figs. 28 and 29). An identical optimum for the sodium periodate concentration was found at 1.8 mM.

20 Similar results were obtained for *S. livingstone* (Figs. 30 and 31). An optimum of 1.8 mM *m*-periodate was found here as well.

Example 11

Extraction of lipopolysaccharides from salmonella spp.

WARNING AND SAFETY PRECAUTIONS

5 *Salmonella* spp., which are used in this extraction method to obtain lipopolysaccharides (LPS), belong to risk class 2 microorganisms. Risk class 2 microorganisms are described as follows: they can induce illnesses like diarrhea and fever, but spreading of the disease is improbable and there are effective prophylaxes or treatment possible. The procedure when working with
10 live bacteria should therefore be carried out with precautions like disinfecting equipment and hands, with 70% ethanol. The waste of living *Salmonella* spp. should be destroyed or disposed of by autoclavation or as biohazard waste, respectively.

LPS compounds are highly pyrogenic and can cause fever. To avoid intake,
15 treat aqueous solutions of LPS with care and wear a mask when working with solid material. If any of these compounds enters the bloodstream, immediately seek medical attention.

0. INTRODUCTION

Salmonella is a gram-negative bacterium, and its outer membrane consists of
20 various antigenic structures, including flagella, outer membrane proteins and Lipopolysaccharides (LPS). The molecule of LPS consists of a so-called lipid A part, which is embedded in the leaflet of the outer membrane, a core region and polysaccharide. The core region is composed of two or three heptoses and two or three residues of eight-carbon, negatively charged monosaccharides
25 KDO. The core region links lipid A to the polysaccharides, which is also known as the O-side chain. This O-side chain is highly variable with respect to its

length and composition between strains, but also within a strain influenced by growth conditions of the *Salmonella*. Despite variation, antigenic structures coded in the PS are unique for a certain *Salmonella* serovar. In fact, antigenic structures O3, O4, O6/7, O8, O9, O10 and O12 represent approximately 90% of known *Salmonella* serovars occurring on porcine products, in particular, Dutch abattoirs.

To detect a humoral response to O antigens as an indication of an exposure of farm animals to *Salomonella*, LPS can be used to probe the binding of raised antibodies to these biomolecules. For this purpose, LPS from *S. typhimurium* (O4,O5,O12), *S. enteritidis* (O9,O12), *S. livingstone* (O6/O7), *S. goldcoast* (O6,O8) and *S. meleagridis* (O3,O10) can be extracted.

An in-house extraction is paramount because LPS from only a limited number of *Salmonella* serovars is commercially available. Furthermore, in-house production can secure a continuous availability of LPS types for a successful antibody detection assay. The in-house extraction method described here for this purpose, is based on a protocol described by Staub (0). Trichloroacetic acid (TCA) extracts LPS containing 1-10% protein contamination. This product is suitable for covalent immobilization of LPS to a carboxymethylated dextran layer coated on a gold metal surface of a biosensor chip (see SOP CHEMIE/A22 (Example 12)). This chip immobilized with LPS, in combination with a Biacore optical SPR biosensor system, can be used to trace *Salmonella*-LPS antibodies in sera also known as serology.

25

1. SCOPE AND FIELD OF APPLICATION

This method describes the extraction of LPS from several *Salmonella* serovars with the use of trichloric acetic acid. Extracted LPS is suitable for

modifications to facilitate its immobilization on a carboxymethylated dextran surface.

2. REFERENCES

- 5 Staub, A.M., Methods in Carbohydrate Chemistry, 5, 92 (1965)
SOP Chemie/A22: Immobilisation of Salmonella-derived LPS onto a biosensor chip (Biacore) and detection of serum antibodies reporting a current or past Salmonella infection (Example 12).
SOP Chemie/A23: Optimisation of protein addition to LPS for immobilization
10 and detection of serum antibodies (Example 13).

3. DEFINITIONS

c = concentration in % (m/v), % (v/v), mol/l or mmol/l as indicated.

15 4. LIST OF ABBREVIATIONS

- 4.1 BGA: Brilliant Green agar
4.2 BHI: Brain Heart Infusion
4.3 BHIa: Brain Heart Infusion agar
4.4 LPS: Lipopolysaccharides
20 4.5 NB: Nutrient broth
4.6 NaCl: Sodium chloride
4.7 NaOH: Sodium hydroxide
4.8 TCA: Trichloric acetic acid

5. PRINCIPLE

Lipopolysaccharides (LPS) are produced by the extraction of *Salmonella* with the aid of trichloroacetic acid (TCA). *Salmonella* is cultured on and then collected from Brain Heart Infusion agar plates. After several washings steps
5 with a saline solution and several centrifugation steps, TCA is added. The acidified suspension is incubated at a low temperature for three hours to solubilise LPS from bacterial cells. The suspension is centrifuged to remove cellular material and the pH is neutralized. LPS is then partly purified and concentrated by ethanol precipitation at low temperature. Finally, salts and
10 ethanol are removed by dialysis, and remaining particles in the retained LPS-containing solution are spun down by centrifugation. The supernatant is lyophilized and weighed to determine the recovery of produced LPS.

6. REAGENTS AND MATERIALS

15 During the procedure, unless stated otherwise, use only reagents of recognized analytical grade and only distilled water or water of equivalent purity. Reference to a company is for information and identification only and does not imply a recommendation unless so stated.

20 6.1 Chemicals

- 6.1.1 Brilliant Green Agar (Oxoid, Basingstroke, England, CM329)
- 6.1.2 Brain Heart Infusion (Oxoid, CM225)
- 6.1.3 Brain Heart Infusion Agar (Oxoid, CM375)
- 6.1.4 Ethanol, absolute (Merck, Darmstad, Germany, 1.00983.2500)
- 25 6.1.5 Glycerol 87% (Merck, 1.04091.1000)
- 6.1.6 Nutrient Broth No 2 (Oxoid, 67)
- 6.1.7 Sodium chloride (Merck, 1.06404.1000)
- 6.1.8 Sodium hydroxide (Merck, 1.06498.1000)

6.1.9 Ethylen glycol (Merck, 9621.2500))

6.1.10 Trichloroacetic acid (Merck, 1.00807.0250)

6.1.11 water was obtained from the Milli Q purification system (8.24)

5 **6.2 *Salmonella* agglutination sera**

6.2.1 anti-O4 (Pro-Lab diagnostics, *Salmonella* reference section of the Central Veterinary Laboratory, Weybridge, Great Britain)

6.2.2 anti-O5 (Pro-Lab diagnostics)

6.2.3 anti-O6,7 (Pro-Lab diagnostics)

10 6.2.4 anti-O8 (Pro-Lab diagnostics)

6.2.5 anti-O9 (Pro-Lab diagnostics)

6.2.6 anti-O12 (Pro-Lab diagnostics)

6.2.7 anti-O Poly A-S (antisera to groups A through S) (Pro-Lab diagnostics)

15 6.2.8 anti-O Poly E (antisera to factors O3, O10, O15, O19, O 34) (Pro-Lab diagnostics)

6.3 Bacterial strains

6.3.1 *Salmonella enteritidis* (#23, phage type 1 strain RIVM, The Netherlands; 90-16-706)

20 6.3.2 *Salmonella goldcoast* (Division's working bank, Utrecht University, The Netherlands)

6.3.3 *Salmonella livingstone* (Division's working bank)

6.3.4 *Salmonella melaegridis* (Division's working bank)

25 6.3.5 *Salmonella typhimurium* X-193 (ASG, Lelystad)

6.4 Reagents

Unless stated otherwise, all prepared media can be stored for 3 months at 2-8°C.

After dispensation of agar media in petri dishes, the plates are cooled at room temperature. After hardening of agar, the plates are inverted and dried at room temperature for two to three days on the tabletop.

- 5 6.4.1 Brilliant Green agar (BGA) plates: Suspend 52 g of BGA (6.1.1) in 1.0 l water (6.1.11). Boil to dissolve the medium completely. Mix well and dispense 15 ml portions in petri dishes (7.29).
- 6.4.2 Brain Heart Infusion (BHI) broth: Dissolve 37 g of BHI broth (6.1.2) in 1.0 l water (6.1.11). Mix well, distribute into final containers and
- 10 sterilize by autoclaving at 121°C for 15 min.
- 6.4.3 Brain Heart Infusion Agar (BHIA): Suspend 47 g BHI agar (6.1.3) in 1.0 l water (6.1.11). Boil to dissolve the medium completely. Mix well and dispense 15 ml portions in petri dishes (7.29).
- 6.4.4 Cooling solution: Mix 1.0 l ethylene glycol (6.1.9) with 3 l water (7.1.11)
- 15 6.4.5 Cold ethanol: Store 1 l ethanol (6.1.4) o/n in a freezer (-18°C) (7.17)
- 6.4.6 Glycerol 87% sterile: Autoclave 50 ml glycerol (6.1.5) at 121°C for 15 min.
- 6.4.7 Nutrient broth (NB): Dissolve 25 g NB (6.1.6) in 1.0 l Water (6.1.11). Mix well, distribute into 100 ml flasks and sterilize by autoclaving at
- 20 121°C for 15 min.
- 6.4.8 Saline ($c = 0.9\%$ (m/v)): Dissolve 9 g sodium chloride (NaCl) (6.1.7) in 1.0 l water (6.1.11). Before use, cool the saline overnight in a refrigerator.
- 6.4.9 Trichloroacetic acid (TCA) solution, $c = 0.25$ mol/l: Dissolve 40.8 g TCA
- 25 (6.1.10) in 1.0 l water (6.1.11).
- 6.4.10 Trichloroacetic acid (TCA) solution, $c = 0.50$ mol/l: Dissolve 81.7 g TCA (6.1.10) in 1.0 l water (6.1.11).
- 6.4.11 Sodium hydroxide (NaOH) solution, $c = 5.0$ mol/l: Dissolve 40 g NaOH (6.1.8) in 200 ml water (6.1.11).
- 30 6.4.12 Sodium hydroxide (NaOH) solution, $c = 0.10$ mol/l: Dilute 1.0 ml $c =$

5 mol/l NaOH (6.4.11) in 49.0 ml water (6.1.11).

7. EQUIPMENT AND APPARATUS

Reference to a company is for information and identification only and does not
5 imply a recommendation unless so stated. Equivalent equipment and
apparatuses may be appropriate as well.

- 7.1 Analytical balance (type AE240, Mettler, Zürich, Switzerland),
weighing range 0 - 40 g (SOP CHEMIE/011)
- 10 7.2 Balance (type MC1 LC2000P, Sartorius Instrumenten BV, Nieuwegein,
The Netherlands)
- 7.3 High speed centrifuge (T-124 Centrikon, Kontron, Zürich, Switzerland)
- 7.4 Superspeed centrifuge (RC-5B, Sorvall, Newton, Connecticut, USA)
- 7.5 Centrifugal polypropylene tubes, for Centrikon 290 ml (253483, Beun
15 de Ronde, Abcoude, The Netherlands)
- 7.6 Centrifugal tubes for rotor SM-24, Sorvall (Sorvall)
- 7.7 Conical tubes (50 ml) (Cellstar 210261, Greiner bio-one, Alphen a/d
Rijn, The Netherlands)
- 7.8 Conical tubes (15 ml) (Cellstar 188271, Greiner bio-one)
- 20 7.9 Cryogenic tubes (1.8 ml) (Cryo's 122261, Greiner bio-one)
- 7.10 Dialysis Membrane clamps (Part#CB-1050, Cellu Sep, Membrane
Filtration Products Inc., Seguin, Texas, USA)
- 7.11 Spectra/Por dialysis tubing MWCO 1000 (Spectrum Laboratories Inc.,
Rancho Dominguez, California, USA)
- 25 7.12 Finn pipette, 200-1000 µl (Labsystems Oy, Helsinki, Finland)
- 7.13 Finn pipette, 40-200 µl (Labsystems Oy, Helsinki, Finland)
- 7.14 Finn pipette, 5-40 µl (Labsystems Oy, Helsinki, Finland)
- 7.15 Freeze dry apparatus (7570001, Labconco Corporation, Kansas City,
Missouri, USA)

- 7.16 Freezer -18°C, freeze range between -16°C to -25°C (Elbanton, Kerkdriel, The Netherlands)
- 7.17 Freezer -80°C (U57085 Classic, New Brunswick Scientific BV, Nijmegen, The Netherlands)
- 5 7.18 Glass beaker, 15 l (Schott-Duran, Omnilabo International BV, Breda, The Netherlands)
- 7.19 Glass slides (AA00000102E, Menzel-Glaser, Braunschweig, Germany)
- 7.20 Ice salt bath: mix approximately 150 g NaCl (7.1.7) with 3 l ice crunch after addition of 150 ml of tap water
- 10 7.21 Home made -4°C incubator (measured temperature of contact liquid to ethanol precipitating LPS extract (9.3.37), approximately 1°C)
 - 7.21.1 Circulation chiller; set temperature -4°C (WK230, Lauda, Lauda-Königshofen, Germany)
 - 7.21.2 Central cooling core PROTEAN II (165-1806, BioRad laboratories BV, 15 Veenendaal, The Netherlands)
 - 7.21.3 Polystyrene box & fitting plastic inner bowl
 - 7.21.3.1 Whole system (8.22.1, 8.22.2) was filled with cooling solution (7.4.4)
- 7.22 Magnetic stirrer (EM3300T, LaboTech B.V., Ochten, The Netherlands)
- 7.23 Milli Q water purification system (Millipore, Bedford, Ma, USA)
- 20 7.24 pH meter (type pH 537, WTW, Weilheim, Germany)
- 7.25 Refrigerator (Elbanton)
- 7.26 Rotor Centrikon (A6.14, Kontron)
- 7.27 Rotor Sorvall (SM-24, Sorvall)
- 7.28 Spatula, Drigalski glass (408014, Omnilabo International BV)
- 25 7.29 Petri dishes (633102, Greiner bio-one)
- 7.30 Powerpette plus (Jencons, Bedfordshire, England)
- 7.31 Vortex-mixer (KS-1, Janke & Kunkel, Staufen, Germany)

8. PROCEDURE

8.1 Preparation stock culture

- 8.1.1 Make an isolate of Salmonella (6.3) by spreading one colony, or the content of an inoculation loop onto a BGA plate (6.4.1).
- 5 8.1.2 Incubate the plate (8.1.1) overnight at 37°C.
- 8.1.3 A single colony is picked from the plate (8.1.2) with an inoculation loop and suspended in 100 ml NB (6.4.7)
- 8.1.4 Incubate overnight at 37°C
- 8.1.5 Add 50 ml glycerol (6.4.5) to 100 ml cultured NB medium (8.1.4)
- 10 8.1.6 Aliquot culture/glycerol mixture (8.1.5) into eleven 14 ml (7.9) (volume of 12.5 ml) and twelve 1.5 ml (7.10) (volume of 1 ml) sterile tubes.
- 8.1.7 One of the 1.5 ml tubes is marked as standard bank. The contents of this tube will only be used to prepare more aliquots of 1 and 12.5 ml via the method here described (8.1).
- 15 8.1.8 Quickly freeze the aliquots at -80°C (8.1.6) and store until use.

8.2 Determination of Salmonella isolate purity by agglutination

- 8.2.1 Pipette 25 µl (7.15) sterile saline (6.4.8) onto a glass slide (7.20).
- 8.2.2 Suspend one colony from the plate cultured Salmonella (8.1.2) in the saline (6.4.8).
- 20 8.2.3 Add a single drop of agglutination serum (6.2) using the facilitated container drop system
- 8.2.4 Mix sera and suspension by gently tilting the slide back and forth for two min.
- 25 8.2.5 Determine agglutination by looking for aggregation formation in front of a black background.
- 8.2.6 Disclose the identity of Salmonella strain by compare the results of aggregation with Table 37.
- 8.2.7 When the identity of the cultured strain differs from predicted

aggregation in Table 37, it can be concluded that the tested culture was not (exclusively) composed of the expected *Salmonella* serotype. In such case, a new stock culture has to be prepared.

8.3 Extraction method

- 5 8.3.1 Prepare 120 BHI agar plates (6.4.3).
- 8.3.2 Thaw a 12.5 ml tube with *Salmonella* stock culture (8.1.8).
- 8.3.3 Spread 100 µl (7.14) stock culture (8.1.8) on each plate (8.3.1) with a spatula (7.29).
- 8.3.4 Incubate overnight at 37°C.
- 10 8.3.5 Fill a 15 l container (7.19) with 10 l water (6.1.11), and cool to 4-8°C in a refrigerator (7.26) until further use.
- 8.3.6 Weigh (7.3) six empty centrifugation tubes (7.6) and note their weight on the quality sheet (see Figure 32.).
- 8.3.7 Put milliQ (6.1.11) and TCA (6.4.9, 6.4.10) containers in prepared ice
15 salt bath (7.21).
- 8.3.8 Take twenty plates (8.3.4) overnight incubated plates and add 1 ml of saline to each plate.
- 8.3.9 Harvest the bacteria by lightly scraping (making round movements) a Drigalski spatula (7.29) over the agar, making a suspension of the
20 bacteria in saline.
- 8.3.10 Collect the suspension in one of the six pre-weighed centrifugation tubes (8.3.6).
- 8.3.11 Wash each plate (8.3.9) twice with 2 ml saline (6.4.8) solution.
- 8.3.12 Combine the suspensions (8.3.11) with the contents of the
25 centrifugation tube (8.3.10).
- 8.3.13 Repeat steps 8.3.8 through 8.3.12 5 times for the remaining plates (8.3.4).
- 8.3.14 Add 100 ml saline (6.4.8) to each centrifugation tube.
- 8.3.15 Tare the tubes (8.3.14)
- 30 8.3.16 Centrifuge for 15 min at 10,000xg and 4°C (7.4, 7.27).

- 8.3.17 Decant the supernatant in a waste container.
- 8.3.18 Suspend each bacterial pellet in 10 ml saline (6.4.8) until a smooth suspension is formed.
- 8.3.19 Add 75 ml saline (6.4.8) to each centrifuge tube.
- 5 8.3.20 Repeat steps 8.3.14 to 8.3.19 once.
- 8.3.21 Repeat steps 8.3.15 to 8.3.17 once.
- 8.3.22 Weigh (7.3) the tubes with bacterial pellet (8.3.21) and note weigh results on the quality sheet (see Figure 32)
- 10 8.3.23 Determine the mass (=m) of the deposited bacteria by subtracting the weight of the empty (8.3.6) with that of the cell-containing tubes (8.3.22).
- 8.3.24 Suspend the bacterial pellets (8.3.22) with x ml (for determination x, see Table 38) of water (6.1.11) by repeatedly drawing in and washing with a 10 ml pipette (7.31).
- 15 8.3.25 Combine the suspensions of two centrifugation tubes in one tube.
- 8.3.26 Repeat 8.3.25 for the remaining centrifugation tubes.
- 8.3.27 Add x ml of y M TCA (6.4.10) (see Table 38 for values x and y).
- 8.3.28 Insert stirring rods in each of the suspensions (8.3.27).
- 8.3.29 Stir the suspensions (8.3.28) for 3 h at 4°C on a magnetic stirrer (7.23).
- 20 8.3.30 Remove the stirring rods.
- 8.3.31 Centrifuge the suspensions for 30 min at 20,000xg and 4°C (7.5, 7.28).
- 8.3.32 Collect the supernatants of the centrifugation tubes in a 500-ml beaker.
- 8.3.33 Adjust the pH of the supernatant with 5 M NaOH (6.4.11) and 0.10 M NaOH (6.4.12) to pH 6.5.
- 25 8.3.34 Determine the volume (v) of the pH adjusted supernatant (8.3.33)
- 8.3.35 Cool the supernatant to freezing point by putting the filled flasks (8.3.34) in a -18°C freezer (7.17) for 30 min
- 8.3.36 Add 2*e ml (for values of e see Table 38) freeze-cold ethanol (6.4.5).
- 30 8.3.37 Cool the solution overnight at -4°C (7.21).

- 8.3.38 Dispense the solution in six centrifuge tubes (7.6).
- 8.3.39 Centrifuge the solution for 30 min at 20,000xg and -4°C (7.4, 7.27).
- 8.3.40 Cut three parts (each 10 cm of length) from the dialysis tube (7.12)
- 8.3.41 Wash the dialysis tubes (8.3.40) briefly with water (6.1.11) and keep
5 them wet in water (6.1.11) until further use
- 8.3.42 Clip a membrane clamp (7.11) at one end of the dialysis tube (8.3.41),
leaving 1 cm tubing free.
- 8.3.43 Decant the supernatant (8.3.39) in a waste bottle.
- 8.3.44 Remove the remaining supernatant from the centrifuge tubes with the
10 help of a pipette (7.13).
- 8.3.45 Add one ml of water (6.1.11) to each centrifuge tube (8.3.44).
- 8.3.46 Suspend the pellets by drawing in and washing out with a one-ml
pipette (7.13).
- 8.3.47 Fill one of the prepared dialysis bags (8.3.42) with the re-suspended
15 pellets (8.3.46) of two centrifuge tubes.
- 8.3.48 Wash each of the tubes with $(z-1)/6$ ml (for the value of z , see Table 38)
water (6.1.11) and combine the wash with the contents of the dialysis
bag (8.3.47).
- 8.3.49 Clip another clamp (6.10) on top of the filled dialysis tube (8.3.48)
20 leaving a small air bubble between solution and clamp.
- 8.3.50 Repeat the steps 8.3.42 to 8.3.49 for the remaining centrifuge tubes.
- 8.3.51 Place the three filled dialysis tubes in pre-cooled water (8.3.5) at 4°C.
- 8.3.52 Incubate dialyse the contents of the tubes (8.3.51) for two days under
continuous gentle stirring conditions at 4°C on a magnetic stirrer
25 (7.23).
- 8.3.53 Refresh the dialysate at least twice, by exchange the water with 7 l
fresh, precooled water (6.1.11).
- 8.3.54 Collect the contents of the dialysis tubes in a 50-ml container (7.8)
- 8.3.55 Divide the collected volume (8.3.54) in centrifuge tubes (7.7)
- 30 8.3.56 Centrifuge the tubes for 20,000xg for 30 min at 4°C (7.5)

- 8.3.57 Weigh an empty 50-ml container (7.8) on an analytic balance (7.2) (without cap) and note its weight on the quality sheet (see Figure 32)
- 8.3.58 Collect the supernatant (8.3.56) in pre-weighed container (8.3.57).
- 8.3.59 Freeze the supernatants in an -80°C freezer (6.15).
- 5 8.3.60 Lyophilize (6.15) the frozen supernatants (8.3.59) until a dry white crystal structure is observed.

NOTE: Handle lyophilized LPS product carefully, as its electrostatic characteristics may easily make it air-borne.

10

NOTE: LPS is a well-known toxin, which when inhaled or swallowed may induce clinical effects in humans. Safety measures, such as masks, should be strictly followed to prevent any intoxication.

- 15 Weigh lyophilized LPS-holding container (8.3.60) on an analytical balance (7.2) and note the resulting mass on the quality sheet (see Figure 32)
- 8.3.62 Calculate the yield of LPS: $(\text{weight LPS (8.3.61)} - \text{weight container (8.3.57)}) / \text{total mass wet cells (8.3.23)} * 100\%$
- 8.3.63 Store lyophilized LPS powder in a closed container at 4°C until further
- 20 use.

An overview of the procedure is depicted in Figure 33.

Table 37. Check board for agglutination readings to disclose the identity of the present bacteria in terms of *Salmonella* species and *Salmonella* serovar

	Agglutination sera							
Salmonella strain	O4 (7.2.1)	O5 (7.2.2)	O6,7 (7.2.3)	O8 (7.2.4)	O9 (7.2.5)	O12 (7.2.6)	poly A-S (7.2.7)	poly E (7.2.8)
<i>S. enteritidis</i> (6.3.1)	-	-	-	-	+	+	+	-
<i>S. goldcoast</i> (6.3.2)	-	-	+	+	-	-	+	-
<i>S. livingstone</i> (6.3.3)	-	-	+	-	-	-	+	-
<i>S. melaegridis</i> (6.3.4)	-	-	-	-	-	-	+	+
<i>S. typhimurium</i> (7.3.4)	+	+	-	-	-	+	+	-

Legend: + = Aggregation formation, agglutination positive

- = No aggregation formation, agglutination negative

Table 38. Conversion table to determine volumes of water (x), TCA (x), ethanol (e) and TCA concentration (y) for extraction and precipitation, and volume of water (z) for dialysis procedures to isolate and purify LPS from *Salmonella*.

Letters noted: m refers to mass (9.3.23), whereas v refers to pH adjusted supernatant volume (9.3.34).

	Extraction		Precipitation	Dialysis
Strain	x ml water (6.1.11) and x ml TCA(6.4.9,6.4.10)	y M TCA	e ml ethanol (6.4.5)	z ml water (6.1.11)
<i>S. enteritidis</i>	$x = m \cdot 5$	0.25 (6.4.9)	$e = 2 \cdot v$	$z = x \cdot 0.1$
<i>S. goldcoast</i>		0.5 (6.4.10)		
<i>S. livingstone</i>		0.5 (6.4.10)		
<i>S. melaegridis</i>		n.e.y.		
<i>S. typhimurium</i>		0.5 (6.4.10)		

n.e.y.: not established yet.

Example 12

Immobilisation of salmonella-derived lps onto a biosensor chip (biacore) and detection of serum antibodies reporting a current or past salmonella infection

5 0. WARNING AND SAFETY PRECAUTIONS

Lipid polysaccharides are potent immunogens, which can bring sensitive persons into a septic shock upon intake or inhalation. Precautions should be made to prevent contact with this biochemical.

10 Sodium periodate is an oxidizing agent and may cause explosions when brought in contact with strong reducing agents.

The procedure utilizes sodium cyanoborohydride. The procedure should therefore be carried out with precautions as with for example hand gloves and a mask. Use this chemical substance only in a chemical fume hood. The material is very toxic to aquatic organisms and may cause long-term adverse
15 effect in the aquatic environment. The material and solution waste should be disposed of as hazardous waste.

1. INTRODUCTION

20 Among the components of the antigenic structure of the genus *Salmonella*, the somatic antigens are important as an instrument to trace immune response in animals upon an invasive infection of this organism. Somatic antigens are located on the polysaccharide part of lipid polysaccharide (LPS), which is a constituent of the bacterial cell wall. After extraction and isolation of carefully chosen LPS (see SOP CHEMIE/A21(Example 11)), antigen-containing LPS is

coupled covalently to a biosensor chip surface to monitor serologically samples for the presence of anti-*Salmonella* antibodies through their binding to the immobilized antigen-containing LPS on the chip surface (see SOP CHEMIE/A23 (Example 13)). This SOP describes the method for oxidation,
5 immobilization of LPS onto the biosensor chip (BIAcore) and the analysis of antibodies in sera.

2. SCOPE AND FIELD OF APPLICATION

To analyze serum samples from chicken for the presence of anti-*Salmonella*
10 antibodies reacting with O4, 5, 6, 7, 8, 9 and 12 somatic antigens.

3. REFERENCES

- Concentration Analysis Handbook, Version AA, December 2001, Biacore AB, Uppsala, Sweden
- 15 BIAapplications Handbook, version AB (reprinted 1998), Biacore
http://www.jp.amershambiosciences.com/tech_support/manual/pdf/dnapuri/52207400af.pdf
- SOP Chemie/A21: Extraction of lipopolysaccharides from *Salmonella* spp. (Example 11)
- 20 SOP Chemie/A23: Optimisation of protein addition to LPS for immobilization and detection of serum antibodies (Example 13).

4. DEFINITIONS

c = concentration in % (m/v), % (v/v), mol/l or mmol/l as indicated.

5. LIST OF ABBREVIATIONS

5.1.1. EDC: *N*-ethyl-*N'*-(3-dimethyl aminopropyl)- carbodiimide hydrochloride

5.1.2. NHS: *N*-hydroxysuccinimide

5 5.1.3. LPS Se: from the *Salmonella* bacteria serovars *enteriditis*

5.1.4. LPS Sg: from the *Salmonella* bacteria serovars *goldcoast*

5.1.5. LPS Sl: from the *Salmonella* bacteria serovars *livingstone*

5.1.6. LPS Sm: from the *Salmonella* bacteria serovars *meleagridis*

5.1.7. LPS St: from the *Salmonella* bacteria serovars *typhimurium*

10 6. PRINCIPLE

LPS is oxidized in the presence of a protein facilitated by sodium periodate. The LPS-protein solution is desalted using a NAP-5 column. The LPS-protein complex is immobilized on a CM5-chip after activation of the carboxymethyl dextran layer on a biosensor chip with the aid of 1-ethyl-3-(3-
15 dimethylaminopropyl)carbodiimide hydrochloride (EDC) and *N*-hydroxysuccinimide (NHS) and carbohydrazide. Bound LPS is then stabilized with sodium cyanoborohydride. Prior to routinely use, the performance of biosensor chip-conjugated LPS to bind anti-*Salmonella* antibodies is assessed using a panel of reference polyclonal agglutination sera.

20 7. REACTIONS

7.1. OXIDATION OF CARBOHYDRATE MOIETY

Periodate will induce an oxidative disruption of linkages between vicinal diols on, in particular, carbohydrate moieties, as in *e.g.* mannose, to yield aldehyde functionalities. This reaction is typically performed in buffers at a pH range
25 between 4.5 and 5.5 in the dark using a freshly prepared 10-100 mM sodium meta-periodate in 0.1 M sodium acetate (Figure 18).

NOTE 1. The positions of conjugations indicated in this scheme, to link up this monosaccharide within a polysaccharide as e.g. LPS are just examples. R' and R indicate the distal and the proximal positions, respectively, in the carbohydrate chain. The oxidation into aldehydes may repeat itself within the polysaccharide chain in each monosaccharide constituent containing vicinal diols.

NOTE 2. Periodate will also oxidize, when present, certain β -aminoethanol derivatives such as the hydroxylysine residues in collagen, as well as methionine (to its sulfoxide) and certain thiols (usually to disulfides). In addition, N-terminal serine and threonine residues of peptides and proteins can be selectively oxidized by periodate to aldehyde groups. These reactions, however, usually occur at a slower rate than oxidation of vicinal diols.

7.2. CONJUGATION TO PROTEIN

Oxidation is performed in the presence of a protein. The bis-aldehyde compounds, like the oxidised monosaccharide constituents in the polysaccharide chain of LPS here above, may react with any amino group in a proteins and may form a Schiff-base linkage resulting in a substituted imine (Figure 19).

NOTE. The substituted imine is stabilised while the complex is attached to the dextran surface, by a reduction facilitated by the cyanoborohydride (reductive animation).

7.3. IMMOBILIZATION TO SENSOR SURFACE

The sensorchip-located carboxymethylated dextran layer is activated by *N*-ethyl-*N*-(3-dimethyl aminopropyl)- carbodiimide hydrochloride (EDC) and *N*-hydroxysuccinimide (NHS). The activation is followed by preparation with carbohydrazine. The reactive aldehyde functionalities react spontaneously with the hydrazide to hydrazones, which are then reduced to stabilise the covalent bonds (Figure 34).

NOTE. The protein R" carries multiple $-NH_2$ groups and can therefore be conjugated with multiple oxidized LPS entities. At the other hand, the polysaccharide part in LPS may carry multiple free aldehyde groups in a single molecule. These aldehyde groups may for a part or complete captured
5 by the hydrazide-dextran layer. The net result may be a very stable complex network of protein-LPS-dextran on the surface.

8. REAGENTS AND MATERIALS

In the complete procedure only reagents of recognized analytical grade and only distilled water or water of equivalent purity are used, unless stated
10 otherwise. Reference to a company is for information and identification only and does not imply a recommendation unless so stated.

8.1. CHEMICALS

- 8.1.1. Acetic acid (J.T. Baker, Deventer, The Netherlands)
- 8.1.2. Amine coupling kit (Biacore AB, Uppsala, Sweden) consisting of:
 - 15 8.1.2.1. Vial containing 115 mg *N*-hydroxysuccinimide (NHS)
 - 8.1.2.2. Vial containing 750 mg 1-ethyl-3-(3-dimethylaminopropyl)carbodiimide hydrochloride (EDC)
 - 8.1.2.3. Vial containing 10.5 ml, $c = 1 \text{ mol/l}$, ethanolamine hydrochloride – sodium hydroxide pH 8.5
- 20 8.1.3. CHAPS (Plus one, Pharmacia Biotech, Uppsala, Sweden)
- 8.1.4. Carbohydrazide, CN_4H_6O (Fluka Chemie GmbH, Buchs, Switzerland)
- 8.1.5. Carboxymethyl-dextran sodium salt (Fluka)
- 8.1.6. Glycine, $c = 10 \text{ mmol/l}$ pH 1.5 (Biacore)
- 25 8.1.7. Guanidine hydrochloride (Calbiochem, San Diego, CA, USA)
- 8.1.8. HBS-EP buffer (Biacore) containing HEPES buffer, $c = 10 \text{ mmol/l}$, pH 7.4, sodium hydrochloride, $c = 150 \text{ mmol/l}$, EDTA, $c = 3 \text{ mmol/l}$ and surfactant P20, $c = 0.005\%(v/v)$.
- 8.1.9. *Salmonella anti group* specific, monoclonal test reagents:

- 8.1.9.1. anti- *Salmonella* gr. B (SIFIN, Berlin, Germany), contains mAb Anti-O4, O5, O27
- 8.1.9.2. anti- *Salmonella* gr. C (SIFIN), contains mAb Anti-O7, O8
- 8.1.9.3. anti- *Salmonella* gr. D (SIFIN), contains mAb Anti-O9, Vi
- 5 8.1.9.4. anti- *Salmonella* gr. E (SIFIN), contains mAb Anti-O3, O19
- 8.1.10. *Salmonella* monovalent 'O' somatic anti sera:
- 8.1.10.1. anti-O4 (Pro-Lab diagnostics, *Salmonella* reference section of the Central Veterinary Laboratory, Weybridge, Great Britain)
- 10 8.1.10.2. anti-O5 (Pro-Lab diagnostics)
- 8.1.10.3. anti-O6,7 (Pro-Lab diagnostics)
- 8.1.10.4. anti-O8 (Pro-Lab diagnostics)
- 8.1.10.5. anti-O9 (Pro-Lab diagnostics)
- 8.1.10.6. anti-O10 (Pro-Lab diagnostics)
- 15 8.1.10.7. anti-O12 (Pro-Lab diagnostics)
- 8.1.10.8. anti-O19 (Pro-Lab diagnostics)
- 8.1.10.9. anti-O Poly E (O3,O10,O15,O19,O34; Pro-Lab diagnostics)
- 8.1.11. *Salmonella* polyvalent 'O' somatic anti sera:
- 8.1.11.1. anti-O Poly A-S
- 20 (O2,O3,O4,O5,O6,7,O8,O9,O10,O11,O12,O13,O15, O16,O17,O18,O19,O20,O21,O22,O23,O28,O30,O34,O35,O38,O 40,O41; Pro-Lab diagnostics)
- 8.1.12. *Salmonella* LPS, in-house isolated LPS by TCA extraction (SOP CHEMIE/A21 (Example 11)) prepared from the *Salmonella*
- 25 bacteria serovars *enteritidis* (Se), *goldcoast* (Sg), *livingstone* (Sl), *meleagridis* (Sm) and *typhimurium* (St) with protein (SOP CHEMIE/A23, example 13)
- 8.1.12.1.1. Aliquots of 0.5 mg LPS are stored at + 4 °C or at lower temperature.
- 30 8.1.13. Avian reference sera

- 8.1.13.1. SPF-CH, SPF serum referred to as negative control serum
(Animal Health Service Ltd. (GD), Deventer, The Netherlands)
- 8.1.13.2. EIA-SE, *Salmonella enteritidis* positive control serum from
chicken for use in ELISA (GD)
- 5 8.1.13.3. EIA-ST, *Salmonella typhimurium* positive control serum
from chicken for use in ELISA (GD)
- 8.1.13.4. SPA-PG, *Salmonella pullorum* positive control serum from
chicken for use in ELISA (GD)
- 8.1.13.5. CH-SI, *Salmonella infantis* positive control serum from
chicken for use in ELISA (GD)
- 10
- 8.1.14. Sodium acetate trihydrate (J.T. Baker, Phillipsburgh, NJ,
USA)
- 8.1.15. Sodium chloride (Merck, Darmstadt, Germany)
- 8.1.16. Sodium cyanoborohydride (NaCNBH_3) (Fluka)
- 15 8.1.17. Sodium hydroxide, $c = 50 \text{ mmol/l}$ (Biacore)
- 8.1.18. Sodium periodate. (Sigma Chemical Comp., St. Louis, MO,
USA)
- NOTE: *This chemical is light-sensitive. Store this material protected
from light*
- 20
- 8.1.19. Triton X-100 (Sigma)
- 8.1.20. Tween 20 (Sigma)
- 8.1.21. Tween 80 (Sigma)
- 25 8.1.22. Water is obtained from a Milli Q water purification system
(MilliQplus)

8.2. SOLUTIONS

- 8.2.1. Acetic acid solution, $c = 0.1 \text{ g/ml}$: Dilute 1 ml acetic acid (8.2) with
9 ml water.

- 8.2.2. Acetate buffer solution, $c = 10$ mmol/l, pH 4.0: dissolve 0.272 g sodium acetate trihydrate (8.2.14) in 180 ml water and adjust to pH 4.0 with acetic acid (8.3.1) and make up to 200 ml with water. This buffer is stable for approximately six months.
- 5 8.2.3. Acetate buffer solution, $c = 1.0$ mol/l, pH 5.5: dissolve 13.6 g sodium acetate trihydrate (8.2.14) in 90 ml water and adjust to pH 5.5 with acetic acid (8.3.1) and make up to 100 ml with water. This buffer is stable for approximately six months.
- 10 8.2.4. Acetate buffer solution, $c = 100$ mmol/l, pH 5.5: dilute 1.0 ml acetate buffer solution, $c = 1.0$ mol/l (8.3.3) with 9.0 ml water.
- 8.2.5. Carbohydrazide solution, $c = 100$ mmol/l: dissolve 9.0 mg carbohydrazide (8.2.4) in 1000 μ l water.
- 8.2.6. Carbohydrazide solution, $c = 5$ mmol/l: dilute 10 μ l carbohydrazide solution (8.3.5) with 190 μ l water. Prepare just before use.
- 15 8.2.7. CHAPS solution, $c = 0.05\%$ (m/v): Dissolve 0.02 g (8.2.3) in 40 ml HBS-EP (8.2.8).
- 8.2.8. Detergents solution, $c = 0.3\%$ (m/v): Dissolve 0.3 g of CHAPS (8.2.3), 0.3 g Tween 20 (8.2.21), 0.3 g Tween 80 (8.2.22) and 0.3 g Triton X-100 (8.2.20) in 100 ml water.
- 20 8.2.9. EDC-solution: reconstitute EDC (8.1.2.2) in 10.0 ml water.
- 8.2.9.1. Fractions of 100 μ l of this solution (8.3.9) are stored in polypropylene tube (9.12) at -18 °C or at lower temperature. The aliquots are stable for two months.
- 25 8.2.9.2. Before use: Thaw frozen aliquots and agitate them gently to ensure homogeneous solutions.
- 8.2.10. Ethanolamine solution: Pipette 200 μ l $c = 1$ mol/l ethanolamine solution (8.1.2.3) in a polypropylene tube (9.12)
- 8.2.11. Guanidine solution, $c = 6$ mol/l: Dissolve 17.18 g guanidine hydrochloride (8.2.7) in 10 ml detergents solution (8.3.8) and adjust
- 30

volume to 30 ml with glycine buffer solution (8.2.6) This solution is stable for approximately three months.

8.2.12. NHS-solution: reconstitute NHS (8.1.2.1) in 10.0 ml water.

5 8.2.12.1. Fractionate 100- μ l aliquots of this solution (8.3.12) in polypropylene tube (9.12). Store at -18°C or at lower temperature. The aliquots are stable for two months.

8.2.12.2. Before use: Thaw frozen aliquots and agitate them gently to ensure homogeneous solutions.

10 8.2.13. Sample dilution buffer: Dissolve 2 g carboxymethyl-dextran sodium salt (8.2.5), 9.97 g sodium chloride (8.2.15) and 0.1 g Tween 80 (8.2.22) in 200 ml HBS-EP (8.2.8).

8.2.14. Sodium cyanoborohydride, $c = 1.00\text{ mol/l}$: dissolve 62.8 mg sodium cyanoborohydride (8.2.16) in 1000 μ l acetate solution pH 4.0 (8.3.2).

15 8.2.15. Sodium cyanoborohydride, $c = 100\text{ mmol/l}$: dilute 20 μ l sodium cyanoborohydride solution (8.3.14) with 180 μ l acetate solution, $c = 10\text{ mmol/l}$, pH 4.0 (8.3.2). Prepare just before use.

20 8.2.16. Sodium hydroxide, $c = 5\text{ mmol/l}$: Dilute 400 μ l sodium hydroxide $c = 50\text{ mmol/l}$ (8.2.17) with 3.6 ml water in glass vial (9.10).

8.2.17. Sodium periodate solution, $c = 100\text{ mmol/l}$: Add to 214 mg sodium periodate (8.2.18) 10.0 ml water.

25 8.2.18. Sodium periodate 'ready to use': Pipette 100 μ l of sodium periodate solution, $c = 100\text{ mmol/l}$ (8.3.17) in a 1.4 ml polypropylene tube (9.13) and dry with a centrifugal evaporator (8.4).

8.2.19. Sodium periodate solution, $c = 50\text{ mmol/l}$: Dissolve sodium periodate 'ready to use' (8.3.18) in 200 μ l acetate solution, $c = 100\text{ mmol/l}$ pH 5.5 (8.3.4). Prepare just before use.

8.3. STANDARD REFERENCE SOLUTION

8.3.1. *Salmonella* anti-O sera: Dilute 20 µl of each serum (8.2.10 and 8.2.11) in 380 µl sample dilution buffer (8.3.13) in a micro titerplate (9.18) with the exception of anti-O5 serum: 2 µl of this serum
5 (8.2.10.2) is diluted in 400 µl sample dilution buffer (8.3.13).

8.3.2. *Salmonella* anti-group specific test reagents: Dilute 4 µl of each serum (8.2.9.1, 8.2.9.2, 8.2.9.3 and 8.2.9.4) in 395 µl sample dilution buffer (8.3.13)

8.3.3. Avian reference sera: Dilute 6 µl sera (8.2.13.1 and 8.2.13.3) in
10 295 µl sample dilution buffer (8.3.13) in a micro titerplate (9.18) and dilute 3 µl sera (8.2.13.2 and 8.2.13.4) in 295 µl sample dilution buffer (8.3.13)

8.3.4. Shake (9.11). Prepare just before use.

15 8.4. AUXILIARY MATERIALS

8.4.1. NAP-5 column (0.5 ml, Sephadex G-25, Amersham Biosciences).

8.4.2. CM5 chips (Biacore).

9. EQUIPMENT AND APPARATUS

9.1. Reference to a company is for information and identification only and
20 does not imply a recommendation unless so stated. Equivalent equipment and apparatuses may be as appropriate as well.

9.2. Analytical balance (type AE 240, Mettler, Zürich, Switzerland)

9.3. Biacore 3000 (Biacore)

9.4. Centrifugal evaporator (Jouan, Saint-Herblain, France)

25 9.5. Finn pipette, 5-40 µl (Lab systems Oy, Helsinki, Finland)

9.6. Finn pipette, 40-200 µl (Lab systems Oy)

9.7. Finn pipette, 200-1000 µl (Lab systems Oy)

9.8. Finn pipette, 1- 5 ml (Lab systems Oy)

- 9.9. Glass collection tube, 5 ml, stoppered (Renes, Zeist, The Netherlands)
- 9.10. Glass vial diameter 16 mm (Biacore)
- 9.11. *Mini shaker with microtiter cup head* (MS 1, Janke & Kunkel, Staufen, Germany)
- 5 9.12. Polypropylene tube, diameter 7 mm, with snap caps (Biacore)
- 9.13. Polypropylene tube, 1.4 ml (Micronic, Lelystad, The Netherlands)
- 9.14. pH meter (type pH 537, WTW, Weilheim, Germany).
- 10 9.15. Sonification bath, Branson 2200 (Branson Ultrasonics B.V., Soest, The Netherlands)
- 9.16. Syringe Filter 0.45 μ m diameter 25 mm (Gelman Sciences, Ann Arbor, MI, USA)
- 15 9.17. Vacuum manifold, to run several NAP-5 columns simultaneously (type SPE-12G, with PTFE stopcock(s), J.T.Baker).
- 9.18. V-bottomed microtiter polystyrene plate, 96-well format (Greiner Bio-one GmbH, Frickenhausen, Germany)
- 9.19. Vortex-mixer (KS-1, Janke & Kunkel)

20 10. SOFTWARE

The biosensor apparatus is operated with Biacore 3000 control software 4.1 (1999-2003).

11. PROCEDURE

11.1. OXIDATION AND DESALTING OF LPS SOLUTION

25

SAFETY PRECAUTION:

Lipid polysaccharides are potent immunogens, which can bring sensitive persons into a septic shock upon intake or inhalation. Precautions should be made to prevent contact with this biochemical.

Sodium periodate is an oxidizing agent and may cause explosions when
5 brought in contact with strong reducing agents.

The procedure utilizes sodium cyanoborohydride. The procedure should therefore be carried out with precautions, such as using hand gloves and a mask. Use this substance only in a chemical fume hood. The material is very toxic to aquatic organisms and may cause long-term adverse effect in the
10 aquatic environment. The material and solution waste should be disposed of as hazardous waste.

11.1.1. OXIDATION

- 11.1.1.1. Add 500 µl acetate buffer pH 5.5 (8.3.4) to the LPS (8.2.12)
15 (See safety precaution)
- 11.1.1.2. Vortex (9.19) thoroughly until the pellet is solved.
- 11.1.1.3. Sonicate (9.15) the solution for 10 minutes and judge the solution for its clearance.
- 11.1.1.4. When clearance is not satisfactory continue sonication until
20 a clear (convalescent) solution is obtained.
- 11.1.1.5. Add 20 µl periodate solution (8.3.19) to the LPS solution (11.1.2.4)
- 11.1.1.6. Vortex (9.19) the solution (11.1.2.5)
- 11.1.1.7. Incubate on ice for 40 min protected from light.
- 25 11.1.1.8. Quench oxidation by desalting the solution (11.1.2.6) as described in 11.1.3

11.1.2. DESALTING

11.1.2.1. Place NAP-5 column(s) (8.5.1) on manifold (9.17).

11.1.2.2. Condition the column(s) (11.1.3.1) by passing three 3-ml portions of acetate buffer (8.3.2) over the column bed on a flow generated by gravity only. Allow the buffer to enter the gel bed completely.

11.1.2.3. Pipette 0.5 ml oxidized LPS solution (11.1.2.8) on the column. Allow the sample to enter the gel bed completely.

11.1.2.4. Elute oxidized LPS with 1 ml of acetate buffer (8.3.2).

Collect eluate in a 5-ml glass tube (9.9).

11.1.2.5. Vortex (9.19) the solution (11.1.3.4) for 10 s.

11.1.2.6. When not immediately used (11.1.3.5) store samples at 4°C to 7°C.

11.1.2.7. Prior to immobilization, the LPS-containing solution is diluted as indicated in the following Table 39.

Table 39.

LPS spp (8.2.12)	µl stock solution (11.1.3.6)	End volume (µl) make with acetate buffer, pH 4.0 (8.3.2)
Se	25	200
Sg	100	200
Sl	100	200
St	100	200
Sm	12.5	200

11.2. IMMOBILIZATION

NOTE: Ensure that all solvent and reagent reservoirs contain sufficient volume to run the method completely before initiating the assay set.

NOTE: Words written in *italic* refer to software commands and menu's etc operating Biacore 3000.

11.2.1. Preparation

- 11.2.1.1. Thaw a portion of EDC (8.3.9.1)
- 5 11.2.1.2. Thaw a portion of NHS solution (8.3.12.1).
- 11.2.1.3. Place the rack Thermo A in the right rack position (R2) and
the Reagent rack in the middle (RR)
- 11.2.1.4. *Command:* Dock a CM5 chip (8.5.2) and prime with HBS-
EP buffer (8.2.8).
- 10 11.2.1.5. Place EDC solution (11.2.1.1) in position R2A1
- 11.2.1.6. Place the NHS solution (11.2.1.2) in position R2A2.
- 11.2.1.7. Place an empty tube (9.12) in position R2A3.
- 11.2.1.8. Place the carbonyldiimidazole solution (8.2.6) in position R2A4.
- 11.2.1.9. Place the ethanolamine solution (8.2.10) in position R2A5.
- 15 11.2.1.10. Place the solution with oxidized LPS (11.1.2.7)
in position R2A6.
- 11.2.1.11. Place the cyanoborohydride solution (8.2.15) in
position R2A7.
- 11.2.1.12. Place the 6 M guanidine solution (8.2.11) in
20 glass vial (9.10) in position RR2
- 11.2.1.13. Place the CHAPS solution (8.2.7) in glass vial
(9.10) in position RR4

PRECAUTION: Place a disposable tube under the waste outlet of the
25 instrument to collect used cyanoborohydride solution. The solution waste
should be disposed of as hazardous waste.

11.2.2. Immobilization of the CM5 chip

11.2.3. *File* (see Figure 35.) → *New application wizard*

11.2.3.1. *Open Template* → Search for file: Wizard immobilisation
(see Figure 36)

5 11.2.3.2. Fill in: Notebook (see Figure 38)

11.2.3.3. *Run, next en start*

11.2.3.4. Note: To see which instruction are in the wizard do *Edit* in
stead off *Run*

10 11.2.4. *Save* the sensorgram. The result files are saved with the
default extension '*.blr*'

11.2.5. The chip is ready for testing *Salmonella* antibodies in sera.

15 **Note:** The immobilization levels obtained should approach the levels as
indicated in Tabel 40. Otherwise consider the immobilization due to oxidation
failure

Table 40. Typical immobilization levels of lps

LPS spp	Immobilization level in RU (Standard deviation) ¹
Se	2365 (959)
Sg	8609 (1969)
Sl	10953 (2135)
St	4836 (1023)

11.3. DETECTION OF ANTI-SALMONELLA ANTIBODIES

20 **Note:** Ensure that all solvent and reagent reservoirs contain
sufficient volume for complete method run of reference sera before
initiating the assay set.

11.3.1. Make the Biacore operational with an appropriate CM5

Example 13

Optimalisation of protein addition to lps for immobilization and detection of serum antibodies

WARNING AND SAFETY PRECAUTIONS

- 5 Lipopolysaccharides (LPS) are highly pyrogenic and can cause fever. To avoid intake, treat aqueous solutions of LPS with care and wear a mask when working with solid material. If any of these compounds enters the bloodstream, directly seek medical attention.

0. INTRODUCTION

- 10 Recent studies indicate that when protein is added to LPS before oxidation, the immobilization to a carboxymethylated dextrane gold layer is made possible and in some cases is improved. Following, serological responses are also made possible and are improved. The optimum of serological responses depends on the percentage of protein added to LPS.
- 15 This protein effect can be obtained by addition of hemoglobin. Hemoglobin is a naturally occurring protein, which can be found in all warm-blooded vertebrates. Therefore cross-reacting anti hemoglobin antibodies in sera are not expected.

1. SCOPE AND FIELD OF APPLICATION

- 20 This method describes the addition of an amount of hemoglobin to lipopolysaccharides (LPS) produced through trichloric acid extraction (see SOP Chemie/A21: Extraction and isolation of Lipopolysaccharides (Example 11)). The optimal hemoglobin percentage is defined as a reaction mixture giving high immobilization levels in combination with maximum serological reaction
- 25 of positive control sera. Furthermore, the production and storage of LPS

reaction mixtures, ready to be used, for immobilization on sensor chips is described.

2. REFERENCES

SOP Chemie/A21: Extraction and isolation of Lipopolysaccharides (version 3;
5 Example 11)

SOP CHEMIE/A22: Immobilization of salmonella derived LPS onto a biosensor chip (BIACORE) and detection of serum antibodies reporting a current or past salmonella infection (version 3; Example12)

3. LIST OF ABBREVIATIONS

- 10 3.1 LPS: Lipopolysaccharides
- 3.2 SOP: Standard Operating Procedure
- 3.3 NaAc: sodium acetate buffer
- 3.4 Hb: Hemoglobin
- 3.5 mQ: milliQ water

15 4. PRINCIPLE

In SOP Chemie/A21 (Example11), the production of *Salmonella spp.* lipopolysaccharides (LPS) is described. Extracted LPS is used as a ligand in an analytical analysis performed on a Biacore 3000 system to detect anti-Salmonella antibodies in sera derived from pigs and chicken (see SOP
20 CHEMIE/A22 (Example 12)). To improve immobilization of LPS, hemoglobin is added before oxidation. To produce a large stock of material to give reproducible immobilization levels, serological data and method performance, LPS is fortified with hemoglobin, divided in aliquots and dried before storage at 4°C. To immobilize a chip, one of the aliquots is batch-wise oxidized and
25 immobilized.

5. REAGENTS AND MATERIALS

During the procedure, unless otherwise stated, use only reagents of recognized analytical grade and only distilled water or water of equivalent purity. Reference to a company is for information and identification only and does not imply a recommendation unless so stated.

5.1 Chemicals

5.1.1 Acetic acid (J.T. Baker, Deventer, The Netherlands)

5.1.2 Hemoglobin, porcine (Sigma-Aldrich, Zwijndrecht, The Netherlands)

10 5.1.3 MilliQ water (0)

5.1.4 Sodium acetate trihydrate (J.T. Baker, Phillipsburgh, NJ, USA)

5.2 *Salmonella* agglutination sera

15 5.2.1 anti-O4 (Pro-Lab diagnostics, Salmonella reference section of the Central Veterinary Laboratory, Weybridge, United Kingdom)

5.2.2 anti-O5 (Pro-Lab diagnostics)

5.2.3 anti-O6,7 (Pro-Lab diagnostics)

5.2.4 anti-O8 (Pro-Lab diagnostics)

20 5.2.5 anti-O9 (Pro-Lab diagnostics)

5.2.6 anti-O10 (Pro-Lab diagnostics)

5.2.7 anti-O12 (Pro-Lab diagnostics)

5.2.8 anti-O19 (Pro-Lab diagnostics)

5.2.9 anti-O Poly A-S (antisera to groups A through S) (Pro-Lab diagnostics)

25 5.2.10 anti-O Poly E (antisera to factors O3, O10, O15, O19, O34) (Pro-Lab diagnostics)

5.3 Group specific *Salmonella* antisera

5.3.1 Enteroclon anti-Salmonella group B (Sifin, Berlin, Germany)

5.3.2 Enteroclon anti-Salmonella group C (Sifin)

5.3.3 Enteroclon anti-Salmonella group D (Sifin)

5.3.4 **Enteroclon anti-Salmonella group E (Sifin)**

5.4 Avian reference sera

- 5 5.4.1 SPF-CH, specific pathogen free (SPF) negative control serum (Animal Health Service Ltd. (GD), Deventer, The Netherlands)
- 5.4.2 EIA-SE, chicken *Salmonella enteritidis* positive control for use in ELISA (GD)
- 10 5.4.3 EIA-ST, chicken *Salmonella typhimurium* positive control for use in ELISA (GD)
- 5.4.4 SPA-PG, chicken *Salmonella pullorum* positive control for use in ELISA (GD)
- 5.4.5 CH-SI, chicken *Salmonella infantis* positive control for use in ELISA (GD)

15 **5.5 Swine reference sera**

- 5.5.1 Sw-Liv, swine *Salmonella livingstone* positive control serum in ELISA (GD)
- 5.5.2 Sw-Typ, swine *Salmonella typhimurium* positive control serum in ELISA (GD)
- 20 5.5.3 Sw-APP, swine *Actinobacillus Pleuropneumoniae* positive control serum in ELISA (GD)

5.6 Lipopolysaccharides

Lipopolysaccharides (LPS) are extracted, lyophilized and stored as described in SOP Chemie/A21 (Example 11).

25

5.6.1 *Salmonella enteritidis* LPS

5.6.2 *Salmonella goldcoast* LPS

5.6.3 *Salmonella livingstone* LPS

5.6.4 *Salmonella meleagridis* LPS

5.6.5 *Salmonella typhimurium* LPS

5.7 Reagents

- 5.7.1 Acetate buffer solution, $c = 10$ mmol/l, pH 4.0: dissolve 0.272 g sodium acetate trihydrate (5.1.4) in 180 ml MQ (5.1.3) and adjust to pH 4.0 with acetic acid (5.1.1) and make up to 200 ml with mQ (5.1.3). This buffer is stable for approximately six months at 4°C.
- 5.7.2 Acetate buffer solution, $c = 1.0$ mol/l, pH 5.5: dissolve 13.6 g sodium acetate trihydrate (5.1.4) in 90 ml MQ (5.1.3) and adjust to pH 5.5 with acetic acid (5.1.1) and make up to 100 ml with mQ (5.1.3). This buffer is stable for approximately six months at 4°C.
- 5.7.3 Hemoglobin stock solution, 5 mg/ml: Dissolve 5 mg hemoglobin (5.1.2) in 1 ml mQ (5.1.3).

5.8 Auxiliary materials

- 5.8.1 CM5 chips (Biacore AB, Uppsala, Sweden).

6. EQUIPMENT AND APPARATUS

Reference to a company is for information and identification only and does not imply a recommendation unless so stated. Equivalent equipment and apparatuses may be appropriate as well.

- 6.1 Centrifugal evaporator (Jouan, Saint-Herblain, France)
- 6.2 Biacore 3000 (Biacore)
- 6.3 Finn pipette, 40-200 µl (Labystems Oy, Helsinki, Finland)
- 6.4 Glass collection tube, 5 ml, 12x75 mm with stop (Renes, Zeist, The Netherlands)
- 6.5 MilliQ installation (Millipore, Bedford, MA, USA)
- 6.6 Sonification bath, Branson 2200 (Branson Ultrasonics, Soest, The Netherlands)
- 6.7 Vortex-mixer (KS-1, Janke & Kunkel, Staufen, Germany)

7. PROCEDURE

7.1 Production of stock solution of lipopolysaccharides

- 5 7.1.1 Collect the produced LPS (see SOP Chemie/A21 (Example 11)) from the refrigerator and let it acclimatize to room temperature.
- 7.1.2 Retrieve the weight of produced LPS (7.1.1) in the tube from the quality data sheet (see SOP Chemie/A21 (Example 11)).
- 7.1.3 Calculate the volume of mQ to be added to LPS using Formulae 1 (9.1).
- 10 7.1.4 Add the calculated volume of mQ (7.1.3) to the LPS tube (7.1.2) (end-concentration LPS: 5 mg/ml).
- 7.1.5 Vortex thoroughly until all powder is solved.
- 7.1.6 Sonicate (6.6) the solution for 10 min and judge the solution for its clearance.
- 15 7.1.7 When clearance (7.1.6) is not satisfactory continue sonication (6.6) until a clear (convalescent) solution is obtained.

7.2 Addition of hemoglobin

- 7.2.1 Prepare four (one for each flowchannel) LPS solution (7.17) dilutions in sodium acetate buffer with variable hemoglobin contents as described in Table 42 (8) in a glass tube (6.4).
- 20 7.2.2 The choice of relative hemoglobin starting amounts added to each newly prepared LPS extraction batch are given in Table 43 (8).
- 7.2.3 Fill up to a total volume of 500 μ l with mQ (5.1.3) as described in Table 42 (8).

25 7.3 Oxidation and desalting (see SOP Chemie/A22; chapter 10.1 (Example 12))

- 7.3.1 Start oxidation from point 10.1.1.2.

7.4 Immobilization (see SOP Chemie/A22, chapter 10.2 (Example 12))

NOTE: Immobilize oxidized LPS fortified with four different relative amounts of hemoglobin (7.3) on one a CM5 chip (5.8.1) so each flow channel represents a different relative amount.

5 7.5 Detection of anti-*Salmonella* antibodies (see SOP Chemie/A22, chapter 10.3 (Example 12))

7.6 Determination of optimal hemoglobin percentage

- 7.6.1 Calculate the mean and standard deviation of the 5 relative responses of each of the agglutination sera listed in (5.2) and the group specific
10 *Salmonella* anti-sera listed in (5.3) per flow channel.
- 7.6.2 Create a clustered column graph with on the x-axis the names of the agglutination and group specific *Salmonella* anti-sera, and on the y-axis the mean of the relative response units for all the different relative amounts of hemoglobin (7.6.1) (see for an example: Figure 41).
- 15 7.6.3 Calculate the mean and standard deviation of the 4 relative responses of each of the avian reference control sera listed in (5.4) and the swine reference sera listed in (5.5) per flow channel.
- 7.6.4 Create a clustered column graph with on the x-axis the names of the avian control and the swine control sera, and on the y-axis the mean of
20 the relative response units for all the different percentages of hemoglobin (7.6.3) (see for an example: Figure 42).
- 7.6.5 Add Y-error bars to both clustered graphs (7.6.2, 7.6.4) by using the standard deviation values for each x-axis column (see for an example: Figures 41 or 42).
- 25 7.6.6 Copy the calculated means of the selected positive expected sera (see Tables 44 to 46 (8)) and the negative SPF chicken sera of all measured relative hemoglobin amounts in a new table (see for an example: Table 47 (8)).
- 7.6.7 Subtract the responses of SPF chicken sera (7.6.6) from the responses

- of the expected positive sera (7.6.6) (see for an example: Table 48 (8)).
- 7.6.8 Determine the highest response per positive serum per relative hemoglobin amount and give this a value of 10 (see for an example Table 49 (8)).
- 5 7.6.9 Calculate for the rest of the flow channels the relative values by using formulae 2a (9.2) (see for an example Table 49 (8)).
- 7.6.10 Calculate the sum of all relative values per hemoglobin percentage (see for an example Table 49 (8)).
- 7.6.11 The optimal hemoglobin percentage (for the four percentages compared) is determined by the highest sum score in the four flow
10 channel/hemoglobin percentages.
- 7.6.12 When the optimal relative hemoglobin amounts (7.6.11) is the highest or lowest hemoglobin amounts compared, steps 7.2 to 7.6.11 have to be repeated with the following conditions.
- 15 7.6.12.1 In case of the lowest hemoglobin amounts for Sg, Sm and Sg (20%) is the most optimal, the amounts of 20% and 30% are repeated in addition to 0% and 10% hemoglobin.
- 7.6.12.2 In case of the highest relative hemoglobin amounts for Se and St is most optimal, the amounts 20 and 30% are repeated in addition to 40
20 and 50% hemoglobin.
- 7.6.12.3 In case the highest relative hemoglobin amounts for Sg, Sm and Sl is most optimal, the percentages 40 and 50% are repeated in addition to 60 and 70% hemoglobin.
- 7.6.13 An optimum of hemoglobin percentage is reached:
- 25 7.6.13.1 When the sum of relative values (7.6.8), per four compared relative hemoglobin amounts, has the highest value.
- 7.6.13.2 In the range of compared hemoglobin where the highest value is detected, a higher and a lower level of hemoglobin addition is also determined.

7.7 Preparation of hemoglobin added vacuum, dried LPS stock

- 7.7.1 The volume of the remaining 5 mg/ml LPS solution after determination of optimal hemoglobin amount is calculated by subtracting tube weight plus LPS solution (7.1.7) by the initial empty tube weight read from the quality data sheet (see SOP Chemie/A21 (Example 11)).
- 7.7.2 The amount of remaining LPS is calculated using Formulae 3 (9.3).
- 7.7.3 Calculate the amount of hemoglobin to be added to LPS using Formulae 4 (9.4).
- 7.7.4 Add the calculated amount of hemoglobin (7.7.1) to the remaining LPS solution to reach an end concentration, which was determined in 7.6.13
- 7.7.5 Invert, vortex and/or sonicate the solution (7.7.4) until the hemoglobin is fully solved.
- 7.7.6 Dispense 100 µl in glass tubes (6.4).
- 7.7.7 Dry the dispensed solution (7.7.6) in a rotating vacuum dryer (6.1) (heating point 1, 15 min., total run time: 60 min.)
- 7.7.8 Stopper the tubes and store at 4 – 7 °C until further use.

Typical baseline responses of *S goldcoast* LPS - hemoglobin complexes immobilized CM5 chip are given in Figure 43.

8. Tables

Table 42. Addition of hemoglobin (5.7.3), sodium acetate buffer (5.7.2) and mQ (5.1.3) to LPS (7.2.1) prior to oxidation.

	Percentage hemoglobin									
	0%	10%	20%	30%	40%	50%	60%	70%	80%	90%
LPS (0)	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
NaAc, 1M pH 5.5 (0)	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl
Hb (5 mg/ml) (0)	0 µl	10 µl	20 µl	30 µl	40 µl	50 µl	60 µl	70 µl	80 µl	90 µl
mQ (0)	350 µl	340 µl	330 µl	320 µl	310 µl	300 µl	290 µl	280 µl	270 µl	260 µl

5

Table 43. Relative hemoglobin starting amounts added to LPS.

PS providing Salmonella strain	hemoglobin			
<i>S.enteritidis</i> (0)	0%	10%	20%	30%
<i>S.goldcoast</i> (0)	20%	30%	40%	50%
<i>S.livingstone</i> (0)	20%	30%	40%	50%
<i>S.meleagridis</i> (0)	20%	30%	40%	50%
<i>S.typhimurium</i> (0)	0%	10%	20%	30%

Table 44. Expected results of (diluted) agglutination serum binding to immobilized Salmonella LPS

LPS providing Salmonella strain	Agglutination sera									
	O4 1:20 (5.2.1)	O5 1:20 0 (5.2.2)	O6,7 1:20 (5.2.3)	O8 1:20 (5.2.4)	O9 1:20 (5.2.5)	O10 1:20 (5.2.6)	O12 1:20 (5.2.7)	O19 1:20 (5.2.8)	O poly A-S 1:20 (5.2.9)	O poly E 1:100 (5.2.10)
<i>S. enteritidis</i> (5.6.1)	-	-	-	-	+	-	+	-	+	-
<i>S. goldcoast</i> (5.6.2)	-	-	+	+	-	-	-	-	+	-
<i>S. livingstone</i> (5.6.3)	-	-	+	-	-	-	-	-	+	-
<i>S. melaegridis</i> (5.6.4)	-	-	-	-	-	+	-	-	+	+
<i>S. typhimurium</i> (5.6.5)	+	+	-	-	-	-	+	-	+	-

Legend: + = positive binding of serum to immobilized LPS

- = no binding of serum to immobilized LPS

⊕ = selected sera to be used in optimal hemoglobin addition determination

Table 45. Expected results of (diluted) avian reference serum binding to immobilized Salmonella LPS

LPS providing Salmonella strain	Avian reference sera				
	CH-SPF 1:20 (5.4.1)	EIA-St 1:20 (5.4.3)	EIA-Se 1:200 (5.4.2)	Spg 1:200 (5.4.4)	Si 1:200 (5.4.5)
<i>S. enteritidis</i> (5.6.1)	-	+	+	+	-
<i>S. goldcoast</i> (5.6.2)	-	-	-	-	+
<i>S. livingstone</i> (5.6.3)	-	-	-	-	+
<i>S. melaegridis</i> (5.6.4)	-	-	-	-	-
<i>S. typhimurium</i> (5.6.5)	-	+	+	+	-

Legend: + = positive binding of serum to immobilized LPS

- = no binding of serum to immobilized LPS

+/- = selected sera to be used in optimal hemoglobin addition determination

Table 46. Expected results of (diluted) swine reference serum binding to immobilized *Salmonella* LPS

LPS providing <i>Salmonella</i> strain	Avian reference sera		
	Sw-Liv (1:20) (5.5.1)	Sw-Typ 1:20 (5.5.2)	Sw-APP 1:20 (5.5.3)
<i>S. enteritidis</i> (5.6.1)	-	+	-
<i>S. goldcoast</i> (5.6.2)	+	-	-
<i>S. livingstone</i> (5.6.3)	+	-	-
<i>S. melaegridis</i> (5.6.4)	-	-	-
<i>S. typhimurium</i> (5.6.5)	-	+	-

Legend: + = positive binding of serum to immobilized LPS

- = no binding of serum to immobilized LPS

+/- = selected sera to be used in optimal hemoglobin addition determination

Table 47. Typical serological responses on *Salmonella goldcoast* immobilized LPS with variable hemoglobin additions (data of two prepared CM5 chips)

Hemoglobin (%)	Immobilization levels (RU)	O6,7 1:20 (5.2.3)	O8 1:20 (5.2.4)	O poly A-S 1:20 (5.2.9)	Anti-Salm Group C (1:100) (5.3.2)	CH-SPF 1:20 (5.4.1)	Sw-Liv (1:20) (5.5.1)
30	4596	273,2	257,1	66,8	76,8	-9,0	53,7
40	4837	262,4	238,3	65,0	67,1	-9,1	51,9
50	6112	271,9	237,5	60,9	71,6	-13,1	47,1
60	9764	221,3	177,8	14,2	41,5	-52,0	10,6
10	2177	179,4	140,0	39,0	55,0	0,5	28,0
20	3469	186,5	143,4	37,4	55,1	-5,3	26,5
30	3690	228,3	186,0	50,1	80,1	-4,7	35,1
40	5922	247,5	198,3	49,7	110,3	-11,2	40,6

BLANK UPON FILING

Table 48. Typical table of subtraction of CH-SPF from the selected positive serological responses (data of Table 47) on *Salmonella goldcoast* immobilized LPS with variable hemoglobin additions (data of two prepared CM5 chips).

Hemoglobin (%)	Immobilization levels (RU)	O6,7 1:20 (5.2.3)	O8 1:20 (5.2.4)	O poly A-S 1:20 (5.2.9)	Anti-Salm Group C (1:100) (5.3.2)	Sw-Liv (1:20) (5.5.1)
30	4596	282,2	266,0	75,8	85,7	62,6
40	4837	271,6	247,4	74,1	76,2	61,0
50	6112	284,9	250,6	74,0	84,6	60,2
60	9764	273,3	229,8	66,2	93,5	62,6
10	2177	178,9	139,5	38,5	54,6	27,5
20	3469	191,8	148,7	42,7	60,4	31,8
30	3690	233,0	190,7	54,8	84,8	39,8
40	5922	258,7	209,5	60,9	121,5	51,8

Table 49. Typical table of relative values score of Table 48 (data of two prepared CM5 chips).

Hemoglobin (%)	Immobilization levels (RU)	O6,7 1:20 (5.3.2)	O8 1:20 (5.2.4)	O poly A-S 1:20 (5.2.9)	Anti-Salm Group C (1:100) (5.3.2)	Sw-Liv (1:20) (5.5.1)	sum
30	4596	10	10	10	9	10	49
40	4837	10	9	10	8	10	47
50	6112	10	9	10	9	10	48
60	9764	10	9	9	10	10	47
10	2177	7	7	6	4	5	30
20	3469	7	7	7	5	6	33
30	3690	9	9	9	7	8	42
40	5922	10	10	10	10	10	50

9. Formulas

9.1 Formulae 1

Calculation of volume of mQ

5

$$v = w/5$$

v = volume of mQ (5.1.3) (in ml)

w = weight of LPS (7.1.2) (mg) (see SOP ChemieA21(Example 11))

10

9.2 Formulae 2

Calculation of relative value of positive sera

15

$$Rv = Mv/Mhv * 10$$

Rv = relative value

Mv = mean value (7.6.1)

Mhv = mean highest value (7.6.8)

20 9.3 Formulae 3

Calculation of amount remaining LPS

$$z = w * 0.005$$

25

z = amount of remaining LPS (in mg)

w = weight of remaining LPS solution (7.7.1) (in mg)

9.4 Formulae 4

Calculation of amount of hemoglobin

$$5 \quad h = y * z * 0.01$$

h = mass of hemoglobin (in mg)

y = optimal hemoglobin percentage (%) (7.6.13)

z = amount of remaining LPS (9.3) (mg)

10

Description of figures

Figure 1. Schematic representation of one embodiment of the method according to the invention

5

Figure 2. Schematic outline of the BIA for detection of bacteria using bacteriophages as indicator organisms. Indicator organisms may be cultured overnight or shorter.

10 **Figure 3.** Reactivity of Hb-fortified, oxidized LPS isolated from *S. typhimurium* (batch St2003.2) with agglutination sera in relative arbitrary biosensor responses (RU). The Hb fortification level of LPS during oxidation is depicted in the figure. The expected binding of the agglutination sera is listed in Table 3.

15

Figure 4. ROC curves from Example 2, Experiment 1. TPF: True-positive fraction;
FPF: false-positive fraction.

20 **Figure 5.** ROC curves from Example 2, Experiment 2. TPF: True-positive fraction;
FPF: false-positive fraction.

Figure 6. Analysis of prepared beads coated with LPS from *S. enteritidis*
25 (reflecting serogroup D), *S. goldcoast* (reflecting serogroup C₂), *S. livingstone* (reflecting serogroup C₁), *S. meleagridis* (reflecting serogroup E) and *S. typhimurium* (reflecting serogroup B). The success of the coating and specificity of the LPS were tested with commercially available monoclonal antisera against O4 (serogroup B), O5 (serogroup B), O7 (serogroup C₁), O8
30 (serogroup C₂) and O9 (serogroup D). Responses are expressed in arbitrary

units as median fluorescence index (MFI) at the Y-axis, whereas the X-axis indicates the type of LPS conjugation of the individual beads.

Figure 7. Analysis of prepared beads coated with LPS reflecting serogroups

5 B, C₁, C₂ and D. The activity of the coating was tested with monoclonal antisera against O4 (serogroup B), O5 (serogroup B), O7 (serogroup C₁), O8 (serogroup C₂) and O9 (serogroup D). Similar to Figure 6, except zoomed in on the lower responses. Notice that response of anti-O5 is under broken. See for details legend of Fig. 6.

10

Figure 8. Comparison of beads coated with two different oxidation batches of oxidized LPS from *S. enteritidis* (reflecting serogroup D) and *S. goldcoast* (reflecting serogroup C₂), *S. livingstone* (reflecting serogroup C₁) and *S. typhimurium* (reflecting serogroup B). The coating was tested with

15 commercially available monoclonal antisera against O5 (serogroup B), O7 (serogroup C₁), O8 (serogroup C₂) and O9 (serogroup D).

Figure 9. Analysis of meat drip and serum from chickens. Commercially available antisera were used to spike meat drip and serum. *Drip*, liquid

20 extract collected from muscle tissue from a chicken, which was tested as *Salmonella*-free using standard ISO methods; *Drip+CH-SPF*, drip that was spiked with serum collected from specific pathogen free (SPF) chickens; *CH-SPF*, serum obtained from specific pathogen free (SPF) chickens; *DripSPA-PG*, drip that was spiked with antiserum reactive with *S. pullorum* and
25 *S. gallinarum*; *SPA-PG*, anti-*S. pullorum* and anti-*S. gallinarum* antiserum; *DripCHSi*, drip that was spiked with chicken serum which was serologically positive for a *S. infantis* infection; *CH-Si*, chicken serum serologically positive for *S. infantis*. The X-axis indicate the type of LPS conjugation of the individual beads. See Figures 6, 7 and 8 for more details.

Figure 10. Analysis of swine sera spiked with commercially available anti-*S. typhimurium* (yellow coloured bars) and anti-*S. livingstone* (cyan coloured bars). In addition, beads in buffer solution (blue coloured bars) and negative swine serum (purple coloured bars) were analysed on beads which were coated with LPS representing serogroups B, C₁, C₂ and D.

Figure 11. Binding of bacteriophage FO1 to immobilized LPS from *S. typhimurium*, *S. enteritidis*, *S. goldcoast* and *S. livingstone* on a Biacore SPR biosensor. PFU, plaque forming units.

Figure 12. Binding of bacteriophage FO1 to an SPR biosensor chip coated with *S. typhimurium* LPS following incubation of *S. typhimurium*, *S. enteritidis*, *S. goldcoast*, *S. livingstone* with 1.2×10^9 PFU bacteriophage FO1. Dotted line indicates the cut off value.

Figure 13. Incubation of different food pathogens and spoilage bacteria in the presence of *Salmonella* spp.-specific bacteriophage FO1. During growth the optical density at λ 600 nm as a measure of bacterial growth was monitored. bl+FO1, blank medium devoid of bacteria supplemented with bacteriophages exclusively.

Figure 14. Incubation of different *Salmonella* serovars in the presence of *Salmonella* spp.-specific bacteriophage FO1. See for more details legend of Fig. 13.

Figure 15. Incubation of different food pathogens and spoilage bacteria in the presence of *Salmonella* spp.-specific bacteriophage FO1. The number of plaque forming units (PFU) was determined following an incubation of 5 h. bl+FO1, blank medium devoid of bacteria supplemented with bacteriophages
5 exclusively.

Figure 16. Incubation of different *Salmonella* serovars in the presence of *Salmonella* spp.-specific bacteriophage FO1. See for more details legend of Fig. 13. FO1stock was not incubated.

10

Figure 17. SPR biosensor analysis of bacteriophage FO1 propagated in different *Salmonella* serovars after concentration and dialysis of the viruses. The suspensions were serial diluted and analysed; the final concentrations of concentrated/diluted bacteriophages is indicated at the X-axis.

15

Figure 18. Oxidation of carbohydrate moiety. R' and R indicate the distal and the proximal positions, respectively, in the carbohydrate chain.

Figure 19. Conjugation to a polyamine containing molecule (R''), such as a
20 protein.

Figure 20. Immobilization to fluorescent beads and stabilization of chemical bonds.

Figure 21. Schematic representation of the procedure of LPS coupling to beads and analysis of serum.

Figure 22a. Examples of Campylobacter-infecting bacteriophages.

5

Figure 22b. Examples of Listeria-infecting bacteriophages.

Figure 22c. Examples of Salmonella-infecting bacteriophages.

10 **Figure 23.** Counting chamber and technique.

Figure 24. Total area of a Bürker-Türk counting chamber (A) of 1 mm², in which B represents the area of 1/16th of the total area.

15 **Figure 25.** Effect of periodate concentration on the immobilization of LPS of *S. enteritidis* on a SPR biosensor chip.

20 **Figure 26.** Effect of periodate concentration on the antigenic activity of immobilized LPS of *S. enteritidis* (batch Se2002.1). LPS was immobilized to a biosensor chip and analyzed in an SPR biosensor (Fig. 25). Range tested was 0.2 mM to 1.8 mM sodium periodate. O9, O12, O poly A-S, S.typh, SE, biosensor response from anti-O9 antisera; anti-O12 antisera, polyclonal antibody against serogroups A to S,

chicken serum positive for *S. typhimurium* and chicken serum positive for *S. enteridis*, respectively.

Figure 27. Effect of periodate concentration on the antigenic activity of
5 immobilized LPS of *S. enteritidis* (batch Se2002.1). Range tested was 1.8 mM to 48.6 mM sodium periodate. See for more details Fig. 26.

Figure 28. Effect of periodate concentration on the immobilization of LPS of
S. goldcoast on a SPR biosensor chip.

10

Figure 29. Effect of periodate concentration on the antigenic activity of
immobilized LPS of *S. goldcoast* (batch Sg2002.1). LPS was immobilized to a
biosensor chip and analyzed in an SPR biosensor (Fig. 28). Range tested was
0.2 mM to 5.4 mM sodium periodate. O6,7, O8, O poly A-S, *S. livingstone*, *S.*
15 *infantis*, biosensor response from anti-O6/7 antisera; anti-O8 antisera, polyclonal
antibody against serogroups A to S, porcine serum positive for *S. livingstone* and chicken
serum positive for *S. infantis*, respectively.

Figure 30. Effect of periodate concentration on the immobilization of LPS of
20 *S. livingstone* on a SPR biosensor chip.

Figure 31. Effect of periodate concentration on the antigenic activity of
immobilized LPS of *S. goldcoast* (batch Sg2002.1). LPS was immobilized to a
biosensor chip and analyzed in an SPR biosensor (Fig. 30). Range tested was
25 0.2 mM to 5.4 mM sodium periodate. O6,7, O8, O poly A-S, *S. livingstone*,

S. infantis, biosensor response from anti-O6/7 antisera; anti-O8 antisera, polyclonal antibody against serogroups A to S, *porcine* serum positive for *S. livingstone* and chicken serum positive for *S. infantis*, respectively.

5 **Figure 32.** Schematic presentation of the procedure.

Figure 33. Quality sheet for LPS extraction process.

10 **Figure 34.** Immobilization of LPS to biosensor surface and stabilization of chemical bonds.

Figure 35. Biacore 3000 control software on COM 1.

15 **Figure 36.** Immobilization wizard.

Figure 37. Immobilization test wizard.

Figure 38. Logging.

20 **Figure 39.** Typical sensorgram of immobilization of oxidized LPS. Report points: 1 baseline; 2 activating EDC/NHS; 3 carbonyldiazide; 4 ethanolamine; 5 Immobilization LPS Se.

25 **Figure 40.** Sensorgram of anti *Salmonella* O Poly A-S analysed on a LPS-containing CM5 chip

Figure 41. Typical serological responses of agglutination sera and group specific *Salmonella* anti sera on a *S. goldcoast* LPS - hemoglobin immobilized CM5 biosensor chip.

30

Figure 42. Typical serological responses of avian reference sera (SPF-CH, EIA St, EIA Se, SPA-PG and CH-Si sera) and swine reference sera (SW-sera) on a *S. goldcoast* LPS - hemoglobin immobilized CM5 biosensor chip.

- 5 Figure 43. Typical baseline responses of a *S. goldcoast* LPS - hemoglobin immobilized CM5 biosensor chip.

References

- Barrow, P., 2000. Serological diagnosis of *Salmonella* by ELISA and other tests. In: Wray, C., Wray, A. (Eds.), *Salmonella in Domestic Animals*. CAB International, Oxon, p. 407.
- Barrow, P.A., Desmidt, M., Ducatelle, R., Guittet, M., van der Heijden, H.M., Holt, P.S., Huis in't Velt, J.H., McDonough, P., Nagaraja, K.V., Porter, R.E., Proux, K., Sisak, F., Staak, C., Steinbach, G., Thorns, C.J., Wray, C., van Zijderveld, F., 1996. World Health Organisation-supervised interlaboratory comparison of ELISAs for the serological detection of *Salmonella enterica* serotype enteritidis in chickens. *Epidemiol. Infect.* 117, 69.
- Berends, B.R., van Knapen, F., Mossel, D.A., Burt, S.A., Snijders, J.M., 1998. Impact on human health of *Salmonella* spp. on pork in The Netherlands and the anticipated effects of some currently proposed control strategies. *Int. J. Food Microbiol.* 44, 219.
- de Vries, N., Zwaagstra, K.A., Huis in't Veld, J.H., van Knapen, F., van Zijderveld, F.G., Kusters, J.G., 1998. Production of monoclonal antibodies specific for the i and 1,2 flagellar antigens of *Salmonella typhimurium* and characterization of their respective epitopes. *Appl. Environ. Microbiol.* 64, 5033.
- Edel, W., van Leeuwen, W.J., Hoogenboom Verdegaal, A.M., 1993. The annual incidence of salmonellosis in humans in The Netherlands. *Tijdschr. Diergeneeskd.* 118, 306.
- Fratamico, P.M., Strobaugh, T.P., Medina, M.B., Gehring, A.G., 1997. Real-time detection of *Escherichia coli* O157:H7 using a surface plasmon resonance

biosensor. Book of Abstracts, 214th ACS National Meeting, Las Vegas, NV, September 7 – 11, AGFD.

- Grimont, P., Grimont, F., Bouvet, P., 2000. Taxonomy of the genus *Salmonella*.
5 In: Wray, C., Wray, A. (Eds.), *Salmonella in Domestic Animals*. CAB International, Oxon, p. 1.

- Holt, P., 2000. Host susceptibility, resistance and immunity to *Salmonella* in animals. In: Wray, C., Wray, A. (Eds.), *Salmonella in Domestic Animals*. CAB
10 International, Oxon, p. 73.

- Hoogenboom Verdegaal, A.M., de Jong, J.C., During, M., Hoogenveen, R., Hoekstra, J.A., 1994. Community-based study of the incidence of gastrointestinal diseases in The Netherlands. *Epidemiol. Infect.* 112, 481.

15

Ivnitski, D., Abdel-Hamid, I., Atanasov, P., Wilkins, E., 1999. Biosensors for detection of pathogenic bacteria. *Biosens. Bioelectron.* 14, 599.

- Jongorius-Gortemaker, B.G.M., B.G., Goverde, R.L., van Knapen, F.,
20 Bergwerff, A.A., 2002. Surface plasmon resonance (BIAcore) detection of serum antibodies against *Salmonella enteritidis* and *Salmonella typhimurium*, *J. Immunol. Meth.* 266, 33–44

- Kamerling JP, Vliegthart JFG, Carbohydrates. In *Clinical Biochemistry, Principles, Methods, Applications. Mass Spectrometry*, edited by Lawson AM
25 (Walter de Gruyter, Berlin, 1989), vol. 1, pp. 175–263.

Kretschmann, E., Raether, H., 1968. Radiative decay of nonradiative surface plasmons excited by light. *Z. Naturforsch. A* 23, 2135.

30

- Liedberg, B., Nylander, C., Lundstroem, I., 1983. Surface plasmon resonance for gas detection and biosensing. *Sens. Actuators* 4, 299.
- Medina, M.B., 1997. SPR biosensor: food science applications. *Food Test. Anal.* 3, 14.
- Medina, M.B., Van Houten, L., Cooke, P.H., Tu, S.I., 1997. Realtime analysis of antibody binding interactions with immobilized *E. coli* O157:H7 cells using the BIAcore. *Biotechnol. Technol.* 11, 173.
- Popoff, M.Y. (2001) *In: Antigenic formulas of the Salmonella serovars*, 8th revision. WHO Collaborating Centre for Reference and Research on *Salmonella*. Institut Pasteur, Paris, France, pp. 150.
- Staub, A.M. (1965) Bacterial lipido-proteino-polysaccharides ('O' somatic antigens): extraction with trichloroacetic acid. In Whistler, R.L. (ed.), *Methods in Carbohydrate Chemistry, Volume V*. Academic Press, New York, pp. 92–93
- Thorns, C.J., Bell, M.M., Sojka, M.G., Nicholas, R.A., 1996. Development and application of enzyme-linked immunosorbent assay for specific detection of *Salmonella enteritidis* infections in chickens based on antibodies to SEF14 fimbrial antigen. *J. Clin. Microbiol.* 34, 792.
- van Asten, A.J., Zwaagstra, K.A., Baay, M.F., Kusters, J.G., Huis in't Veld, J.H., van der Zeijst, B.A., 1995. Identification of the domain which determines the g,m serotype of the flagellin of *Salmonella enteritidis*. *J. Bacteriol.* 177, 1610.
- Van Pelt, W., Van de Giessen, A., Van Leeuwen, W., Wannet, W., Henken, A., Evers, E., De Wit, M., Van Duynhoven, Y., 1999. Oorsprong, omvang en kosten

van humane salmonellose. Deel 1. Oorsprong van humane salmonellose met betrekking tot varken, rund, kip, ei en overige bronnen. Infectieziekten Bull. 10, 240.

- 5 W. van Pelt, W.J.B. Wannet, A.W. van de Giessen, Y.T.H.P. van Duynhoven, 2003. Trends in gastroenteritis in the Netherlands Lowest incidences in 2002 ever; the lull before the storm? Infectieziekten bull. 14, 424.

Wilkons S.G., 1996, Bacterial lipopolysaccharides- Themes and variations.

- 10 Progress in Lipid Research, volume 35, issue 3, sept, p. 283-343

Yamane, Y., Awamura, N., Fujii, H., Ohta, H., Toyota, Y., Otsuki, K., Inoue, T., 2000. Establishment of an enzyme-linked immunosorbent assay with a coated deflagellated Salmonella enteritidis antigen for detection of a specific
15 chicken antibody. Avian Dis. 44, 291.

References Example 2

- Baay, M.F., Huis in 't Veld, J.H. (1993) Alternative antigens reduce cross-
20 reactions in an ELISA for the detection of Salmonella enteritidis in poultry. J. Appl. Bacteriol. 74, 243.

- Bergwerff, A.A., van Knapen, F. (2006) Surface Plasmon Resonance Biosensors for Detection of Pathogenic Micro-organisms: Strategies to Secure Food and
25 Environmental Safety. J. A.O.A.C. Int., accepted for publication.

Bergwerff, A.A., van Knapen, F. (2003) Sensing pathogens: Screening strategies in food and environmental safety. Biacore Journal 2, 10-15.

- Bokkers, E.G.M. (2002). The Action Plan Salmonella Poultry Layer Sector, 2001+. Dutch Product Board for Livestock, Meat and Eggs. Available from http://bedrijfsnet.pve.agro.nl/pls/pbs/docs/folder/BEDRIJFSNET_US_CA/HOOFDTHEMAS/KWALITEIT_VOEDSELVEILIGHEID/SALMONELLA_CAMPYLOBACTER/ACTIEPLANNEN/EIERSECTOR/ARTIKEL_ACTIEPLANSALMONELLA_%20EIERSECTOR2001.PDF. In Dutch. Last visited: December 2005
- Desmidt, M., Ducatelle, R., Haesebrouck, F., de Groot, P.A., Verlinden, M., Wijffels, R., Hinton, M., Bale, J.A., Allen, V.M. (1996) Detection of antibodies to *Salmonella enteritidis* in sera and yolks from experimentally and naturally infected chickens. *Vet. Rec.* 138, 223-6.
- Van Duijkeren E., Wannet, W.J., Houwers, D.J., van Pelt, W. (2002). Serotype and phage type distribution of salmonella strains isolated from humans, cattle, pigs, and chickens in the Netherlands from 1984 to 2001. *J. Clin. Microbiol.* 40, 3980-5.
- Van Duynhoven, Y.T., de Jager, C.M., Kortbeek, L.M., Vennema, H., Koopmans, M.P., van Leusden, F., van der Poel, W.H., van den Broek, M.J. (2005) Explosie Project Team. A one-year intensified study of outbreaks of gastroenteritis in The Netherlands. *Epidemiol. Infect.* 133, 9-21.
- Fischer, I.S.T. (2004) International trends in *Salmonella* serotypes 1998-2003 – a surveillance report from the Enter-net international surveillance network. *Eurosurveillance* 9, 9-10.
- Gast, R.K., Nasir, M.S., Jolley, M.E., Holt, P.S., Stone, H.D. (2002) Detection of experimental *Salmonella enteritidis* and *S. typhimurium* infections in laying hens by fluorescence polarization assay for egg yolk antibodies. *Poult. Sci.* 81, 1128-31.

- Gast, R.K., Porter, R.E. Jr., Holt P.S. (1997) Applying tests for specific yolk antibodies to predict contamination by *Salmonella enteritidis* in eggs from experimentally infected laying hens. *Avian Dis.* 41, 195-202.
- 5
- Gast, R.K., Beard, C.W. (1991) Detection of *Salmonella* serogroup D-specific antibodies in the yolks of eggs laid by hens infected with *Salmonella enteritidis*. *Poult. Sci.* 70, 1273-6.
- 10
- Guerin, M.T., Martin, S.W., Darlington, G.A., Rajic, A. (2005). A temporal study of *Salmonella* serovars in animals in Alberta between 1990 and 2001. *Can. J. Vet. Res.* 69,88-99.
- Hassan, J.O., Barrow, P.A., Mockett, A.P., McLeod, S. (1990) Antibody
- 15 response to experimental *Salmonella typhimurium* infection in chickens measured by ELISA. *Vet. Rec.* 126, 519.
- Holt, P.S., Stone, H.D., Gast, R.K., Greene, C.R. (2000) Application of the agar gel precipitin test to detect antibodies to *Salmonella enterica* serovar
- 20 *enteritidis* in serum and egg yolks from infected hens. *Poult. Sci.* 79, 1246-50.
- Jongerijs-Gortemaker, B.G., Goverde, R.L., van Knapen, F., Bergwerff, A.A. (2002) Surface plasmon resonance (biacore) detection of serum antibodies against *Salmonella enteritidis* and *Salmonella typhimurium*. *J. Immunol. Meth.* 266, 33-44.
- 25
- Li, Z.Z., Gong, F.C., Shen, G.L., Yu, R.Q. (2002) Bacteria-modified amperometric immunosensor for a *Brucella melitensis* antibody assay. *Anal. Sci.* 18, 625-30.

Liu, G.D., Wu, Z.Y., Wang, S.P., Shen, G.L., Yu, R.Q. (2001) Renewable amperometric immunosensor for *Schistosoma japonicum* antibody assay. *Anal. Chem.* 73, 3219-26.

- 5 Metz, C.E., Herman, B.A., Roe, C.A. (1998) Statistical comparison of two ROC curve estimates obtained from partially-paired datasets. *Med. Decis. Making* 18, 110-121

- Van Pelt, W., van de Giessen, A., van Leeuwen, W., Wannet, W., Henken, A.,
10 Evers, E., de Wit, M., van Duynhoven, Y. (1999) Oorsprong, omvang en kosten van humane salmonellose. Deel 1. Oorsprong van humane salmonellose met betrekking tot varken, rund, kip, ei en overige bronnen. *Infectieziekten Bull.* 10, 240-243. In Dutch.

- 15 Proux, K., Jouy, E., Houdayer, C., Protais, J., Dibb-Fuller, M., Boscher, E., Gillard, A., Gracieux, P., Gilbert, F., Beaumont, C., Duchet-Suchaux, M. (2002) Reliable ELISAs showing differences between resistant and susceptible lines in hens orally inoculated with *Salmonella enteritidis*. *Vet. Res.* 33, 23-33.

- 20 Pyrohova, L.V., Starodub, M.F., Artiukh, V.P., Nahaieva, L.I., Dobrosol, H.I. (2002) Express diagnostics of bovine leucosis by immune sensor based on surface plasmon resonance. *Ukr. Biokhim. Zh.* 74, 88-92. In Ukrainian.

- Sachsenweger, O., Lohr, J.E., Kusters, J. (1994) Evaluation of three
25 commercial ELISA test kits for the detection of antibodies against *Salmonella enteritidis*. *Tierarztl. Prax.* 22, 350-7. In German.

- Skov, M.N., Feld, N.C., Carstensen, B., Madsen, M. (2002) The serologic response to *Salmonella enteritidis* and *Salmonella typhimurium* in
30 experimentally infected chickens, followed by an indirect lipopolysaccharide

enzyme-linked immunosorbent assay and bacteriologic examinations through a one-year period. *Avian Dis.* 46, 265-73.

- 5 Su, X., Low, S., Kwang, J., Chew, V.H.T., Li, S.F.Y. (2001). Piezoelectric quartz crystal based veterinary diagnosis for *Salmonella enteritidis* infection in chicken and egg. *Sensors and Actuators B: Chemical* 75, 29-35.

- 10 Sunwoo, H.H., Nakano, T., Dixon, W.T., Sim, J.S. (1996) Immune responses in chickens against lipopolysaccharide of *Escherichia coli* and *Salmonella typhimurium*. *Poult. Sci.* 75, 342-5.

- 15 Uttenthaler, E., Kosslinger, C. and Drost, S. (1998) Characterization of immobilization methods for African swine fever virus protein and antibodies with a piezoelectric immunosensor. *Biosens. Bioelectron.* 13, 1279-86.

Vetcha, S., Wilkins, E., Yates, T. (2002) Detection of hantavirus infection in hemolyzed mouse blood using alkaline phosphatase conjugate. *Biosens. Bioelectron.* 17, 901-9.

- 20 De Vries, N., Zwaagstra, K.A., Huis in't Veld, J.H., van Knapen, F., van Zijderveld, F.G. and Kusters J.G. (1998) Production of monoclonal antibodies specific for the i and 1,2 flagellar antigens of *Salmonella typhimurium* and characterization of their respective epitopes. *Appl. Environ. Microbiol.* 64, 5033.

- 25 Van Zijderveld, F.G., van Zijderveld-van Bommel, A.M. and Anakotta, J. (1992) Comparison of four different enzyme-linked immunosorbent assays for serological diagnosis of *Salmonella enteritidis* infections in experimentally infected chickens. *J. Clin. Microbiol.* 30, 2560-6.

Zweig, N.H. and Campbell, G. (1993) Receiver-operating characteristic (ROC) plots: a fundamental evaluation tool in clinical medicine. Clin. Chem. 39, 561-77. Erratum in: Clin Chem 1993, 39, 1589.

5 References Example 3

Benoit P.W. and D.W. Donahue. 2003. Methods for rapid separation and concentration of bacteria in food that bypass time-consuming cultural enrichment. J. Food Prot. 66(10):1935-1948.

10

Che Y., Y. Li, and M. Slavik. 2001. Detection of *Campylobacter jejuni* in poultry samples using an enzyme-linked immunoassay coupled with an enzyme electrode. Biosens. Bioelectron 16(9-12):791-797.

15 Christensen B., H. Sommer, H. Rosenquist and N. Neilsen. 2001. Risk assessment on *Campylobacter jejuni* in chicken product. Available Source: http://www.foodrisk.org/risk_assessments.cfm?item1=Campylobacter&item2=All+Commodities&Submit=Submit, June 15, 2004.

20 Coker A.O. 2000. Incidence, trends and sources of Campylobacteriosis in developing countries - An overview, pp. 44-48. In Proceedings of a WHO Consultation of Experts, Copenhagen, Denmark, 21-25 November 2000. Available Source: <http://www.who.int/emcdocuments/zoonoses/whocdscsraph20017c.html>, June 15, 2004.

25

FSAI (Food Safety Authority of Ireland). 2002. Control of *Campylobacter* species in the food chain. Available source: http://www.fsai.ie/publications/report/Campylobacter_report.pdf, June 15, 2004.

Lake, R., A. Hudson, P. Cressey and G. Nortje. 2003. Risk Profile: *Campylobacter jejuni/coli* in poultry (whole and pieces). Client Report FW0109. Christchurch: ESR. (Institute of Environmental Science & Research Limited). Available source: <http://www.nzfsa.govt.nz/science-technology/risk-profiles/campylobacter.pdf>, June 15, 2004.

Stern, N.J. and J.E. Line. 2000. *Campylobacter*, pp. 1040-1056. In Lund B.M., T.C. Baird-Parker and G.W. Gould (eds.). The Microbiological Safety and Quality of Food Volume II. Aspen Publishers Inc., Gaithersburg, Maryland.

10

Tauxe R.V. 2000. Major risk factors for human *Campylobacteriosis*—An overview, pp. 65-66. In Proceedings of a WHO Consultation of Experts, Copenhagen, Denmark, 21-25 November 2000. Available Source: <http://www.who.int/emc-documents/zoonoses/whocdscsgraph20017c.html>, June 15, 2004.

15

Waller D.F. and S.A. Ogata. 2000. Quantitative immunocapture PCR assay for detection of *Campylobacter jejuni* in foods. Appl. Environ. Microbiol. 66(9): 4115-4118.

20

Yu L.S., J. Uknalis, and S.I. Tu. 2001. Immunomagnetic separation methods for the isolation of *Campylobacter jejuni* from ground poultry meats. J. Immunol. Methods. 256(1-2): 11-18.

25 References Example 4

Salin pork: Final Report of Salmonella in Pork (SALINPORK) Preharvest and harvest control options based on epidemiology, diagnostic and economic research (2000) EU project FAIR1 CT95-0400, (Lo Fo Wong, D.M.A. & Hald, T, eds), Copenhagen, Denmark, pp. 251.

30

- European Commission - Health & Consumer Protection Directorate-General
(2004) Trends and sources of zoonotic agents in animals, feedingstuffs, food
5 and man in the European Union and Norway in 2002, Part 1, SANCO/29/2004,
pp. 32

References Examples 6, 7 and 8

- 10 Lindberg, A A. Bacterial surface carbohydrates and bacteriophage adsorption.
In: Sutherland I E. , editor; Sutherland I E. , editor. Surface carbohydrates of
the prokaryotic cells. London, United Kingdom: Academic Press; 1977. pp.
289–356.
- 15 Lindberg AA, Holme T. 1969. Influence of O side chains on the attachment of
the Felix O-1 bacteriophage to Salmonella bacteria. J Bacteriol. 99(2):513-9.
- Hirsh, Dwight C. and Martin, Lori D. 1983. Detection of Salmonella spp. in
milk by using Felix-O1 bacteriophage and high-pressure liquid
20 chromatography. Appl Environ. Microbiol. 46(5):1243-5.

2006237730 05 Jul 2011

CLAIMS:

1. A method for immobilisation of a polysaccharide on a detection system carrier, comprising contacting said polysaccharide with an oxidising agent in the presence of a polymer comprising at least two amine and/or amide groups to obtain a polysaccharide-polymer complex
5 and coupling said polysaccharide-polymer complex to said carrier.
2. A method according to claim 1, said polymer is protein.
3. A method according to claim 1 or claim 2, wherein said polysaccharide is derived from a gram-negative bacterium.
4. A method according to any one of claims 1 to 3, wherein said polysaccharide is derived
10 from an enterobacteriaceae.
5. A method according to any one of claims 1 to 4, wherein said polysaccharide is derived from a salmonella (sub)species.
6. A method according to any one of claims 1 to 5, wherein said polysaccharide is a lipopolysaccharide.
- 15 7. A method according to any one of claims 2 to 6, wherein said protein, is haemoglobin or myoglobin.
8. A method according to any one of claims 1 to 7, wherein said oxidising agent is (sodium) *m*-periodate.
9. A method according to any one of claims 1 to 8, further comprising activating the surface
20 of said carrier.
10. A method according to any one of claims 1 to 9, wherein said carrier comprises a glass surface coated with gold.
11. A method according to anyone of claims 1 to 10, wherein said carrier is modified with a coating comprising a carboxyl group donor, preferably a carboxymethylated dextran layer.
- 25 12. A method according to claim 11, wherein said carboxyl group donor and/or dextran layer is activated with 1-ethyl-3-(3-dimethylaminopropyl)carbodiimide hydrochloride, N-hydroxysuccinimide and carbonylhydrazide .

22 Feb 2012

2006237730

13. A method according to any one of claims 1 to 12, comprising at least two different polysaccharides.
14. A method according to anyone of claims 1 to 13, wherein said carrier is a biosensor chip.
15. A detection system carrier obtained by the method according to any one of claims 1 to 14.
- 5 16. A detection system carrier comprising an immobilised polysaccharide-protein complex on its surface, wherein said carrier comprises a coating comprising a carboxyl group donor covalently linked to a polysaccharide, wherein said carboxyl group donor and said polysaccharide are covalently linked to each other via a polymer comprising at least two amine and/or amide groups, wherein at least said polysaccharide is linked to said polymer via a
0 periodate oxidised vincinal diol on said polysaccharide and an amine and/or amide group on said polymer, wherein said polymer is a protein.
17. A carrier according to claim 16, wherein the carboxyl group donor is a carboxymethylated dextran.
18. A carrier according to claim 16 or 17, which is a microsphere or bead.
- 15 19. A carrier according to claim 18, wherein said microsphere or bead is a polystyrene microsphere or bead.
20. A carrier according to any one of claims 15-19, that is coded.
21. A carrier according to claim 20, wherein said carrier is coded by the presence of a label.
22. A carrier according to claim 21, wherein said label comprises a fluorescent or
20 phosphorescent colour.
23. A collection of microsphere or beads comprising at least two differently encoded microsphere or beads according to any one of claims 20-22.
24. A collection of microsphere or beads according to claim 23, wherein each of said differently encoded microsphere or beads comprises a polysaccharide that comprises a different
25 antigen.
25. A biosensor comprising a carrier according to claim 15 or 22.

2006237730 05 Jul 2011

26. A Surface Plasmon Resonance detection system comprising a biosensor according to claim 25.

27. A method for determining the presence of an antibody directed to an antigen of a gram-negative bacteria in a sample, comprising contacting said sample with a carrier according to any one of claims 15 or 22 or a biosensor according to claim 25 and determining whether the carrier has bound any antibody.

28. A method according to claim 27, wherein said sample is blood, blood-derived liquid material, tissue-derived fluids, such as meat drip, milk, egg, fluids from an eye, saliva or faeces.

29. A method for determining the presence of a gram-negative bacterium in a sample, comprising contacting said sample with a predetermined amount of antibodies directed against an antigen of said bacterium and determining the amount of antibodies not bound to said bacterium with a carrier according claim 15 or 22 or a biosensor according to claim 25.

30. A method according to any one of claims 27 to 29, wherein binding to said carrier or said biosensor is determined by Plasmon Surface Resonance.

31. A method according to any one of claims 27 to 30, wherein said sample is obtained from a human or an animal.

32. A method for determining the presence of a gram-negative bacterium in a sample comprising

- contacting said sample with target bacteria-specific, bacteriophages and allowing the bacteriophages to infect said sample

- removing non-bound and/or non-invading bacteriophages resulting in a bacteriophage infected sample

- bringing the bacteriophage infected sample into contact with an indicator organism susceptible for the used bacteriophages

- incubate during at least one bacteriophage multiplication cycle

- recover the bacteriophages to obtain a bacteriophage-containing sample

- analyse said bacteriophage-containing sample with a carrier according to claim 15 or 22 or a biosensor according to claim 25.

33. A method according to claim 32, wherein said sample is obtained from a human, a plant or an animal.

34. A method according to claim 32 and 33, wherein said bacteriophage comprises a bacteriophage of figure 22a, 22b and/or 22c.

35. A carrier according to any one of claims 15-22 comprising a bacteriophage of figure 22a, 22b and/or 22c.

36. A method according to claim 1, substantially as hereinbefore described.

2006237730 05 Jul 2011

1/41




Fig. 1

2/41




Fig. 2




Fig. 3

4/41




Fig. 4

5/41




Fig. 5

6/41




Fig. 6

7/41




Fig. 7

8/41




Fig. 8

9/41




Fig. 9

10/41




Fig. 10

11/41




Fig. 11

12/41

Binding of FO1 after Salmonella incubation determined on LPS-St: 2420 RU




Fig. 12

13/41




Fig. 13

14/41




Fig. 14

15/41




Fig. 15

16/41




Fig. 16

17/41




Fig. 17

18/41




Fig. 18

19/41




Fig. 19




Fig. 20

20/41




Fig. 21

21/41

NTCC12669 , NTCC12670 , NTCC12671 , NTCC12672 , NTCC12673 ,
NTCC12674 , NTCC12675 , NTCC12676 , NTCC12677 , NTCC12678 ,
NTCC12679 , NTCC12680 , NTCC12681 , NTCC12682 , NTCC12683 ,
NTCC12684

Fig. 22a


A005 , A006 , A020 , A500 , A502 , A511 , A118 , A620 , A640 , B012 , B021 ,
B024 , B025 , B035 , B051 , B053 , B054 , B055 , B056 , B101 , B110 , B545 ,
B604 , B653 , C707 , D441 , HSO47 , H1OG , H8/73 , H19 , H21 , H43 , H46 ,
H107 , H108 , H110 , H163/84 , H312 , H340 , H387 , H391/73 , H684/74 ,
H924A , PSA , U153 , fMLUP5 , (syn= P35) , 00241 , 00611 , 02971A , 02971C ,
5/476 , 5/911 , 5/939 , 5/11302 , 5/11605 , 5/11704 , 184 , 575 , 633 , 699/694 ,
744 , 900 , 1090 , 1317 , 1444 , 1652 , 1806 , 1807 , 1921/959 , 1921/11367 ,
1921/11500 , 1921/11566 , 1921/12460 , 1921/12582 , 1967 , 2389 , 2425 , 2671 ,
2685 , 3274 , 3550 , 3551 , 3552 , 4276 , 4277 , 4292 , 4477 , 5337 , 5348/11363 ,
5348/11646 , 5348/12430 , 5348/12434 , 10072 , 11355C , 11711A , 12029 ,
12981 , 13441 , 90666 , 90816 , 93253 , 907515 , 910716

Fig. 22b


b , Beccles , CT , d , Dundee , f , Fels 2 , GI , GIII , GVI , GVIII , k , K , i , j , L ,
O1 , (syn= O-1) , (syn= O1) , (syn= O-I) , (syn= 7) , O2 , O3 , P3 , P9a , P10 ,
Sab3 , Sab5 , San15 , San17 , SI , Taunton , ViI , (syn= ViI) , 9

Fig. 22c

22/41



Pic. 6



Pic. 6.1

Fig. 23

23/41




Fig.24

24/41




Fig. 25

25/41




Fig. 26

26/41




Fig. 27

27/41




Fig. 28

28/41




Fig. 29

29/41




Fig. 30

30/41




Fig. 31

31/41

Quality sheet of LPS extraction						
Strain:	Salmonella			Inoculation date/time		
date:				harvesting date/time		
batch:				time of incubation	hours	
amount of BHIa						
Collecting cells + extraction:						
weight centrifugated culture (g)			TCA extraction			
tube #	mass tube empty	mass tube + pellet	mass cells (m)	volume x milli Q	volume x TCA	
tube						
tube						
tube						
tube						
tube						
tube						
		total mass cells		+		
Determination of x:		$x = m \cdot 5$				
molarity y TCA (M):		y =				
incubation time TCA				hours		
ethanolprecipitation:						
Adjust pH		> 6.5?	yes / no	max. pH		
		start pH				
		end pH				
volume supernatant/erlenmeyer (v)				ml		
volume ethanol absoluut:		$e = 2 \cdot v$		ml		
Dialysis:						
volume milli Q:		$z = 0.1 \cdot x$		ml		
Lyophilization:						
		mass tube empty	mass tube + LPS	mass LPS		
yield LPS:				g		
recovery LPS:		(mass LPS/mass total cells) * 100%			%	
REMARKS:						

Fig. 32

32/41




Fig. 33

33/41




Fig. 34




Fig. 35




Fig. 36

35/41




Fig. 37

36/41

running buffer: HBS-EP lot
 chip: CM5 lot
 immob date:
 reagents:
 EDC lot: Reconstitution date:
 NHS lot: Reconstitution date:
 ethanolamine lot:
 Preparation date carbonyldiimidazole solution:
 Preparation date cyanoborohydride solution:
LPS:
 concentration.: µg/ml
 dilution before immobilization:
 Immobilization on:
 Fc1:
 oxidation date: Immobilization level RU
 Fc2:
 oxidation date: Immobilization level RU
 Fc 3:
 oxidation date: Immobilization level RU
 Fc4:
 oxidation date: Immobilization level RU
 wash:
 3 x 5 µl guanidine solution (**Error! Reference source not found.**) and 3 x 5 µl CHAPS (**Error! Reference source not found.**)
 sample diluted with:
 HBS-EP + 0.85 M NaCl + 1.0% CM dextraan + 0.05% Tween 80
 flow 20 µl/min
 injection volume 40 µl
 Salmonella 'O'somatic anti sera:

antiserum:	diluted	lot	verd.
O4	1:20		20+380
O5	1:200		2+400
O6,7	1:20		20+380
O8	1:20		20+380
O9	1:20		20+380
O10	1:20		20+380
O12	1:20		20+380
O19	1:20		20+380
O poly E	1:20		20+380
O poly A-S	1:20		20+380
antiSalm grB	1:100		4+395
antiSalm grC	1:100		4+395
antiSalm grD	1:100		4+395
antiSalm grE	1:100		4+395
Chicken sera:	1:50 or 1:200		
SPF-CH	1:50		6+295
EIA-ST	1:50		6+295
EIA-SE	1:200		1.5+300
SPA-PG	1:200		1.5+300
CH-SI	1:200		1.5+300

 wash:
 *10 µl guanidine solution (**Error! Reference source not found.**) and
 flow 100 µl/min

Fig. 38

37/41




Fig. 39




Fig. 40




Fig. 41




Fig. 42

41/41




Fig. 43

SUBSTITUTE SHEET (RULE 26)