

(12) STANDARD PATENT
(19) AUSTRALIAN PATENT OFFICE

(11) Application No. **AU 2009240505 B2**

(54) Title
Isoprene synthase variants for improved microbial production of isoprene

(51) International Patent Classification(s)
C12N 5/10 (2006.01) **C12N 15/52** (2006.01)
C12N 9/88 (2006.01)

(21) Application No: **2009240505** (22) Date of Filing: **2009.04.23**

(87) WIPO No: **WO09/132220**

(30) Priority Data

(31) Number	(32) Date	(33) Country
61/125,336	2008.04.23	US

(43) Publication Date: **2009.10.29**

(44) Accepted Journal Date: **2013.09.05**

(71) Applicant(s)
DANISCO US INC.;THE GOODYEAR TIRE & RUBBER COMPANY

(72) Inventor(s)
Cervin, Marguerite A.;Whited, Gregory M.;Miasnikov, Andrei;Peres, Caroline M.;Weyler, Walter;Bott, Richard R.;Wells, Derek H.;Rife, Christopher;McAuliffe, Joseph;Kellis, James

(74) Agent / Attorney
Houlihan2, Level 1 70 Doncaster Road, BALWYN NORTH, VIC, 3104

(56) Related Art
WO 2008/003078 A2, 3 January 2008

(19) World Intellectual Property Organization
International Bureau(43) International Publication Date
29 October 2009 (29.10.2009)

PCT

(10) International Publication Number
WO 2009/132220 A9

(51) International Patent Classification:

C12N 9/88 (2006.01) C12N 5/10 (2006.01)
C12N 15/52 (2006.01)

(21) International Application Number:

PCT/US2009/041581

(22) International Filing Date:

23 April 2009 (23.04.2009)

(25) Filing Language:

English

(26) Publication Language:

English

(30) Priority Data:

61/125,336 23 April 2008 (23.04.2008) US

(71) Applicants (for all designated States except US):
DANISCO US INC. [US/US]; 925 Page Mill Road, Palo Alto, CA 94304-1013 (US). **THE GOODYEAR TIRE & RUBBER COMPANY** [US/US]; 1144 East Market Street, Akron, OH 44316-0001 (US).

(72) Inventors; and


(75) Inventors/Applicants (for US only): **CERVIN, Marguerite, A.** [CA/US]; 404 Palomar Drive, Redwood City, CA 94062 (US). **WHITED, Gregory, M.** [US/US]; 304 South Road, Belmont, CA 94002 (US). **MIASNIKOV,****Andrei** [FI/US]; 1769 Latham Street, #a, Mountain View, CA 94041 (US). **PERES, Caroline, M.** [BE/US]; 2328 Santa Catalina Street, Palo Alto, CA 94303 (US). **WEYLER, Walter** [US/US]; 252 Theresa Street, San Francisco, CA 94112 (US). **WELLS, Derek, H.** [US/US]; 775b Loma Verde Avenue, Palo Alto, CA 94303 (US). **BOTT, Richard, R.** [US/US]; 3032 Hillside Drive, Burlingame, CA 94010 (US). **RIFE, Christopher, Lee.** [US/US]; 1499 Hudson Street, Redwood City, California 94061 (US).(74) Agents: **SHIEH-NEWTON, Terri** et al.; Morrison & Forster LLP, 755 Page Mill Road, Palo Alto, CA 94304-1018 (US).

(81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BR, BW, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IS, JP, KE, KG, KM, KN, KP, KR, KZ, LA, LC, LK, LR, LS, LT, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PG, PH, PL, PT, RO, RS, RU, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.

[Continued on next page]

(54) Title: ISOPRENE SYNTHASE VARIANTS FOR IMPROVED MICROBIAL PRODUCTION OF ISOPRENE

Figure 92



(57) Abstract: The present invention provides methods and compositions comprising at least one isoprene synthase enzyme with improved catalytic activity and/or solubility. In particular, the present invention provides variant plant isoprene synthases for increased isoprene production in microbial host cells. Biosynthetically produced isoprene of the present invention finds use in the manufacture of rubber and elastomers.



(84) Designated States (*unless otherwise indicated, for every kind of regional protection available*): ARIPO (BW, GH, GM, KE, LS, MW, MZ, NA, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, SE, SI, SK, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

— *with sequence listing part of description (Rule 5.2(a))*

(88) Date of publication of the international search report:
21 January 2010

(48) Date of publication of this corrected version:
6 May 2010

(15) Information about Correction:
see Notice of 6 May 2010

Published:

— *with international search report (Art. 21(3))*

ISOPRENE SYNTHASE VARIANTS FOR IMPROVED MICROBIAL PRODUCTION OF ISOPRENE

CROSS-REFERENCE TO RELATED APPLICATIONS

[0001] This application claims the benefit of U.S. provisional application serial no. 61/125,336 filed April 23, 2008, which is hereby incorporated by reference in its entirety.

FIELD OF THE INVENTION

[0002] The present invention provides methods and compositions comprising at least one isoprene synthase enzyme with improved catalytic activity and/or solubility. In particular, the present invention provides variant plant isoprene synthases for increased isoprene production in microbial host cells. Biosynthetically produced isoprene of the present invention finds use in the manufacture of rubber and elastomers.

BACKGROUND OF THE INVENTION

[0003] Isoprenoids are isoprene polymers that find use in pharmaceuticals, nutraceuticals, flavors, fragrances, and rubber products. Natural isoprenoid supplies, however, are limited due to ecological concerns. For this reason, and to provide isoprenoid compositions having fewer impurities and greater uniformity, isoprenoids such as rubber are often produced synthetically.

[0004] Isoprene (2-methyl-1,3-butadiene) is a volatile hydrocarbon that is insoluble in water and soluble in alcohol. Commercially viable quantities of isoprene can be obtained by direct isolation from petroleum C5 cracking fractions or by dehydration of C5 isoalkanes or isoalkenes (Weissermel and Arpe, *Industrial Organic Chemistry*, 4th ed., Wiley-VCH, pp. 117-122, 2003). The C5 skeleton can also be synthesized from smaller subunits. It would be desirable, however, to have a commercially viable method of producing isoprene that was independent of nonrenewable resources.

[0005] Biosynthetic production of isoprene occurs by two distinct metabolic pathways (Julsing et al., *Appl Microbiol Biotechnol*, 75:1377-1384, 2007). In eukaryotes and archae, isoprene is formed via the mevalonate (MVA) pathway, while some eubacteria and higher plants produce isoprene via the methylerythritol phosphate (MEP) pathway. Isoprene emissions from plants are light and temperature-dependent with increases linked to leaf development. An isoprene-producing enzyme, isoprene synthase, has been identified in Aspen trees (Silver and Fall, *Plant Physiol*, 97:1588-1591, 1991; and Silver and Fall, *J Biol*

2009240505 12 Jul 2013

Chem, 270:13010-13016, 1995) and is believed to be responsible for the *in vivo* production of isoprene from whole leaves. Bacterial production of isoprene has also been described (Kuzma et al., Curr Microbiol, 30:97-103, 1995; and Wilkins, Chemosphere, 32:1427-1434, 1996), and varies in amount with the phase of bacterial growth and the nutrient content of the culture medium (U.S. Patent No. 5,849,970 to Fall et al.; and Wagner et al., J Bacteriol, 181:4700-4703, 1999, both herein incorporated by reference in their entirety). The levels of isoprene obtainable through bacterial systems of the prior art, however, are insufficient for commercial uses.

[0006] Thus what the art needs is an efficient, large scale, bacterial isoprene production process to provide feedstock for the manufacture of isoprenoids.

[0007] All patents, patent applications, articles and publications mentioned herein are hereby expressly incorporated herein by reference.

SUMMARY OF THE INVENTION

[0008] The present invention provides methods and compositions comprising at least one isoprene synthase enzyme with improved catalytic activity and/or solubility. In particular, the present invention provides variant plant isoprene synthases for increased isoprene production in microbial host cells. Biosynthetically produced isoprene of the present invention finds use in the manufacture of rubber and elastomers.

[0008A] The present invention provides an isolated isoprene synthase variant comprising one or more amino acid substitution(s) at one or more amino acid positions corresponding to a poplar isoprene synthase having the sequence of SEQ ID NO:120 wherein said substitution(s) are selected from the group consisting of V10M, F12S, T15A, E18G, V58I, V58F, L70Q, L70R, L70V, L70T, T71P, V79L, E89D, G94A, S119F, F120L, G127R, E175V, T212I, S257A, R262G, A266G, F280L, N297K, F305L, L319M, E323K, A328T, D342E, A359T, K366N, E368D, L374M, S396T, V418S, K438N, H440R, T442A, I449V, A469S, K500R, K505Q, G507S, S509N, F511Y, and N532K; and wherein the variant is capable of more effectively converting dimethylallyl diphosphate (DMAPP) to isoprene, as compared to an isoprene synthase variant without a substitution.

[0009] For example, the present invention provides isolated isoprene synthase variants, wherein the variant comprises a substitution at a position corresponding to one or more residues (one, two, three, four, five, six, seven, eight, nine or ten) of a kudzu isoprene synthase comprising the amino acid sequence set forth in SEQ ID NO: 2. In some embodiments, the isoprene synthase variant is a kudzu (*Pueraria sp.*) isoprene synthase

2009240505 12 Jul 2013

variant or a poplar (*Populus sp.*) isoprene synthase variant. In some embodiments, the one or more residues are selected from but not limited to the group consisting of L26, E30, F31, Q33, L35, E36, N37, L39, K40, V41, K43, L44, R61, V62, D63, Q65, K87, E94, N95, L99, D100, N105, K137, E138, G143, E144, N182, L184, K185, G187, N189, T190, P225, H226, K247, T257, E258, M259, D266, N334, D353, S357, I358I, E361, N389, I392, I393, K398, E401, C421, Q423, Q424, E425, D426, H430, L432, R433, S434, D437, R443, L462, E463, H476, N478, D479, Q485, D508, P513, A515, Q532, Y533, L537, G538, R539, Y542, A543, and P557. In some embodiments, the one or more residues are selected from but not limited to the group consisting of P24, N25, Y309, D310, L377, F381, E384, Y399, N402, A403, S406, S407, G409, A411, L413, F449, A456, T457, S458, A459, A460, E461, L462, E463,

R464, G465, E466, T467, T468, N469, M523, S527, and Y531. In some embodiments, the one or more residues are selected from but not limited to the group consisting of A20, N21, Y22, Q23, R271, W278, F299, V302, and S408. The present invention also provides an isolated isoprene synthase variant having an A20G substitution in a kudzu isoprene synthase having the amino acid sequence set forth in SEQ ID NO: 2. In a subset of these embodiments, the variant comprises at least two substitutions (two, three, four, five, six, seven, eight, nine or ten), wherein one of the substitutions is an A20G substitution in a kudzu isoprene synthase having the amino acid sequence set forth in SEQ ID NO: 2. The present invention also provides an isolated isoprene synthase variant having an S408D substitution in a kudzu isoprene synthase having the amino acid sequence set forth in SEQ ID NO: 2. In a subset of these embodiments, the variant comprises at least two substitutions (two, three, four, five, six, seven, eight, nine or ten), wherein one of the substitutions is an S408D substitution in a kudzu isoprene synthase having the amino acid sequence set forth in SEQ ID NO: 2. In some preferred embodiments, the isoprene synthase variant has at least one improved property as compared to wild-type isoprene synthase. In some particularly preferred embodiments, the at least one improved property is selected from but not limited to the group consisting of specific activity (production of isoprene from dimethylallyl diphosphate), and solubility.

[0010] In addition, the present invention further provides a polynucleotide sequence encoding the isoprene synthase variant. Also provided is an expression vector comprising a polynucleotide sequence encoding the isoprene synthase variant in operable combination with a promoter. In further embodiments, the present invention provides a host cell comprising the expression vector. Also provided is a lysate of the host cell, wherein the lysate further comprises lysozyme. In some embodiments, the lysate has a neutral pH (6.5 to 7.5), while in other embodiments the lysate has a basic pH (above 7.5 and below 9.5). The present invention also provides methods of producing isoprene, comprising: (a) providing host cells comprising the expression vector; and (b) culturing the host cells under conditions suitable for producing isoprene. In some embodiments, the methods further comprise (c) recovering the isoprene. In still further embodiments, the methods further comprise (d) polymerizing the isoprene. The present invention further provides methods of detecting isoprene synthase activity, comprising: (a) culturing host cells comprising the expression vector under conditions suitable for producing the isoprene synthase variant; (b) lysing the host cells with a lysis buffer comprising lysozyme to produce a cell lysate; and (c) detecting isoprene synthase activity in the cell lysate by measuring isoprene production from dimethylallyl diphosphate (DMAPP). In some embodiments, the host is selected from but not limited to the

group consisting of gram-positive bacterial cells, gram-negative bacterial cells, filamentous fungal cells, and yeast cells. In some preferred embodiments, the host is selected from but not limited to the group consisting of *Escherichia sp.* (*E. coli*), *Pantoea sp.* (*P. citrea*), *Bacillus sp.* (*B. subtilis*), *Yarrowia sp.* (*Y. lipolytica*), and *Trichoderma* (*T. reesei*). In some embodiments, the host cells are cultured in a medium that includes a carbon source selected from but not limited to the group consisting of glucose, glycerol, glycerine, dihydroxyacetone, yeast extract, biomass, molasses, sucrose, and oil.

[0011] Moreover, the present invention provides methods of detecting isoprene in a plurality of samples (high-throughput screening), comprising: (a) providing: i) a plurality of samples each comprising isoprene synthase; ii) a glass plate comprising a plurality of wells; and iii) a seal for the glass plate; (b) placing the plurality of samples in the plurality of wells of the glass plate; (c) sealing the glass plate with the seal to produce a sealed glass plate having a headspace associated with the sample in each of the plurality of wells; (d) incubating the glass plate under conditions in which the isoprene synthase is active; and (e) detecting isoprene in the headspace. In some embodiments, the isoprene is detected by gas chromatography-mass spectrometry (GC-MS). In some embodiments, the plurality of samples comprise host cells comprising an expression vector comprising a polynucleotide sequence encoding an isoprene synthase variant in operable combination with a promoter. In some embodiments, the plurality of samples comprise a lysate of the host cells, lysozyme, and dimethylallyl diphosphate (DMAPP). In some preferred embodiments, the glass plate is a deep-well glass block. In some preferred embodiments, the plurality of wells comprises at least 24 wells (preferably at least 48 wells, more preferably at least 96 wells, still more preferably at least 192 wells, and most preferably at least 384 wells). In particularly preferred embodiments, the plurality of wells each comprise a volume of 2 ml or less (preferably 2 ml to 0.2 ml).

[0012] Additionally the present invention provides a host cell comprising a heterologous polynucleotide sequence encoding an isoprene synthase variant in operable combination with a promoter, wherein the isoprene synthase variant comprises a substitution at a position corresponding to one or more residues (one, two, three, four, five, six, seven, eight, nine or ten) of a kudzu isoprene synthase comprising the amino acid sequence set forth in SEQ ID NO: 2. In some embodiments, the one or more residues are selected from but not limited to the group consisting of L26, E30, F31, Q33, L35, E36, N37, L39, K40, V41, K43, L44, R61, V62, D63, Q65, K87, E94, N95, L99, D100, N105, K137, E138, G143, E144, N182, L184, K185, G187, N189, T190, P225, H226, K247, T257, E258, M259, D266, N334, D353, S357,

I358I, E361, N389, I392, I393, K398, E401, C421, Q423, Q424, E425, D426, H430, L432, R433, S434, D437, R443, L462, E463, H476, N478, D479, Q485, D508, P513, A515, Q532, Y533, L537, G538, R539, Y542, A543, and P557. In some embodiments, the one or more residues are selected from but not limited to the group consisting of P24, N25, Y309, D310, L377, F381, E384, Y399, N402, A403, S406, S407, G409, A411, L413, F449, A456, T457, S458, A459, A460, E461, L462, E463, R464, G465, E466, T467, T468, N469, M523, S527, and Y531. In some embodiments, the one or more residues are selected from but not limited to the group consisting of A20, N21, Y22, Q23, R271, W278, F299, V302, and S408. The present invention also provides an isolated isoprene synthase variant having an A20G substitution and/or an S408D substitution in a kudzu isoprene synthase having the amino acid sequence set forth in SEQ ID NO: 2. In some preferred embodiments, the isoprene synthase variant has at least one improved property as compared to wild-type isoprene synthase. In some particularly preferred embodiments, the at least one improved property is selected from but not limited to the group consisting of specific activity (production of isoprene from dimethylallyl diphosphate), and solubility. In some preferred embodiments, the polynucleotide sequence is contained within a plasmid. In other preferred embodiments, the polynucleotide sequence is integrated into a chromosome of the host cell. In some embodiments, the host is selected from but not limited to the group consisting of gram-positive bacterial cells, gram-negative bacterial cells, filamentous fungal cells, and yeast cells. In some preferred embodiments, the host is selected from but not limited to the group consisting of *Escherichia sp.* (*E. coli*), *Pantoea sp.* (*P. citrea*), *Bacillus sp.* (*B. subtilis*), *Yarrowia sp.* (*Y. lipolytica*), and *Trichoderma (T. reesei)*. In some embodiments, the host cells are cultured in a medium that includes a carbon source selected from but not limited to the group consisting of glucose, glycerol, glycerine, dihydroxyacetone, yeast extract, biomass, molasses, sucrose, and oil. In some embodiments, the host cell further comprises a heterologous or native nucleic acid encoding an IDI polypeptide and/or a heterologous or native nucleic acid encoding a DXS polypeptide, sometimes in combination with the native DXP pathway (for example, expression of *dxs* and *idi* in *E. coli* in addition to the native DXP pathway). Alternatively the entire DXP pathway (Figure 15) maybe expressed on a plasmid or integrated on the chromosome as an operon, with a single promoter controlling expression, or promoters of varying strengths (example GI 1.20, GI 1.5, or GI 1.6) controlling one or more of the individual genes. In some embodiments, the host cell further comprises one or more nucleic acids encoding an IDI polypeptide and a DXS polypeptide, while in some preferred embodiments, one vector encodes the isoprene synthase variant, the IDI

polypeptide, and the DXS polypeptide. In some embodiments, the host cell further comprises a heterologous nucleic acid encoding an MVA pathway polypeptide (e.g., an MVA pathway polypeptide from *Saccharomyces cerevisia* or *Enterococcus faecalis*). In some embodiments, the host cell further comprises one or more nucleic acids encoding an MVA pathway polypeptide and a DXS polypeptide, while in some preferred embodiments, one vector encodes the isoprene synthase variant, the MVA pathway polypeptide, and the DXS polypeptide. In some preferred embodiments, the host cell further comprises one or more nucleic acids encoding a DXS polypeptide, an IDI polypeptide, or one or more of the rest of the DXP pathway polypeptides, and a MVA pathway polypeptide. In some embodiments, the vector further comprises a selectable marker (e.g., antibiotic resistance nucleic acid). Also provided are methods of producing isoprene, comprising: (a) culturing the host cells under suitable culture conditions for production of isoprene; and (b) producing the isoprene. In some embodiments, the methods further comprise (c) recovering the isoprene. In some preferred embodiments, the methods further comprise (d) polymerizing isoprene. The present invention also provides methods of producing isoprene synthase, comprising: (a) providing: (i) a host cell; and (ii) a nucleic acid encoding an isoprene synthase variant in operable combination with a promoter, wherein the isoprene synthase variant comprises a substitution at a position corresponding to one or more residues (one, two, three, four, five, six, seven, eight, nine or ten) of a kudzu isoprene synthase comprising the amino acid sequence set forth in SEQ ID NO: 2;

(b) contacting the host cell with the nucleic acid to produce a transformed host cell; and (c) culturing the transformed host cells under suitable culture conditions for production of isoprene synthase.

[0013] In another aspect, the invention provides for isolated poplar isoprene synthase variants. In one embodiment, the variant comprises a truncation in the N-terminal portion of isoprene synthase. In another embodiment, the isoprene synthase variant has an increased specific activity compared to a full length isoprene synthase. In another embodiment, the isoprene synthase is *P. alba* isoprene synthase of SEQ ID NO:120. In another embodiment, wherein the variant is selected from the group consisting of: an MEA variant (SEQ ID NO:122), an MSV variant (SEQ ID NO:124), an MVS variant (SEQ ID NO:126), an MTE variant (SEQ ID NO:128), an MNV variant (SEQ ID NO:130). In another embodiment, the variant is an MEA variant (SEQ ID NO:122). In another embodiment, the variant is selected from the group consisting of: a TRC (-3) variant (SEQ ID NO:136), a TRC (-4) variant (SEQ ID NO:138), a TRC (-5) variant (SEQ ID NO:140), a TRC (-6) variant (SEQ ID NO:142) and

a TRC (-7) variant (SEQ ID NO:144). In another embodiment, the variant is a MET variant of *P. tremuloides* isoprene synthase (SEQ ID NO:146). In another embodiment, the variant is a MET variant of *P. trichosarpa* isoprene synthase (SEQ ID NO:148).

[0014] In another aspect, the invention provides for isolated poplar isoprene synthase variants, wherein the variant comprises a substitution of one or more amino acid residues of a wild type isoprene synthase; and wherein the isoprene synthase variant has increased isoprene synthase activity compared to a wild type isoprene synthase. In one embodiment, the increased isoprene synthase activity is indicated by a host cell comprising the isoprene variant growing at a faster rate in the presence of dimethylallyl pyrophosphate (DMAPP) compared to a host cell comprising a parent isoprene synthase. In another embodiment, the isoprene synthase is the *P. alba* isoprene synthase of SEQ ID NO:120. In another embodiment, the variant comprises one or more amino acid substitutions selected from the group consisting of V10M, F12S, T15A, E18G, V58I, V58F, L70Q, L70V, L70T, T71P, V79L, E89D, G94A, S119F, F120L, G127R, E175V, T212I, S257A, R262G, A266G, F280L, N297K, F305L, L319M, E323K, A328T, D342E, A359T, K366N, E368D, L374M, S396T, V418S, K438N, H440R, T442I, T442A, I449V, A469S, K500R, K505Q, G507S, S509N, F511Y, and N532K. In another embodiment, at least one amino acid substitution is a L70R substitution. In another embodiment, the variant comprises one or more amino acid substitutions selected from the group consisting of G127R/F511Y, L70Q/G94A/R262G/F305L, F12S/T15A/E18G/N297K, S396T/T442I, V10M/E323K, F120L/A266G, K438N/K500R, V79L/S509N, E175V/S257A/E368D/A469S, T71P/L374M, F280L/H440R, E89D/H440R, V58F/A328T/N532K, S119F/D342E/I449V, and K366N/G507S.

[0015] In another aspect, the invention provides for a crystalline form of a polypeptide comprising the amino acid residues of SEQ ID NO:120 (Figure 19).

[0016] In another aspect, the invention provides for methods of producing isoprene, comprising: (a) providing a host cell comprising an expression vector comprising a polynucleotide sequence encoding an isoprene synthase variant; and (b) culturing the host cell under conditions suitable for producing isoprene. In one embodiment, the method further comprises (c) recovering the isoprene. In another embodiment, the method further comprises (d) polymerizing the isoprene.

[0017] In another aspect, the invention provides for methods of detecting isoprene synthase activity, comprising: (a) culturing a host cell comprising the expression vector under conditions suitable for producing an isoprene synthase variant; (b) lysing the host cells with a

lysis buffer comprising lysozyme to produce a cell lysate; and (c) detecting isoprene synthase activity in the cell lysate by measuring isoprene production from dimethylallyl diphosphate (DMAPP). In one embodiment, the host cell is selected from the group consisting of gram-positive bacterial cells, gram-negative bacterial cells, filamentous fungal cells, and yeast cells. In another embodiment,

the host cell is selected from the group consisting of *Escherichia sp.* (*E. coli*), *Pantoea sp.* (*P. citrea*), *Bacillus sp.* (*B. subtilis*), *Yarrowia sp.* (*Y. lipolytica*), and *Trichoderma* (*T. reesei*).

In another embodiment, the host cell is cultured in a medium that includes a carbon source selected from the group consisting of glucose, glycerol, glycerine, dihydroxyacetone, yeast extract, biomass, molasses, sucrose, and oil.

[0018] In another aspect, the invention provides for host cells comprising a heterologous polynucleotide sequence encoding an isoprene synthase variant in operable combination with a promoter, wherein the isoprene synthase variant comprises a substitution at a position corresponding to one or more residues (one, two, three, four, five, six, seven, eight, nine or ten) of a poplar isoprene synthase. In one embodiment, the isoprene synthase is the *P. alba* isoprene synthase of SEQ ID NO:120. In another embodiment, the variant is selected from the group consisting of: an MEA variant (SEQ ID NO:122), an MSV variant (SEQ ID NO:124), an MVS variant (SEQ ID NO:126), an MTE variant (SEQ ID NO:128), an MNV variant (SEQ ID NO:130). In another embodiment, the variant is selected from the group consisting of: a TRC (-3) variant (SEQ ID NO:136), a TRC (-4) variant (SEQ ID NO:138), a TRC (-5) variant (SEQ ID NO:140), a TRC (-6) variant (SEQ ID NO:142) and a TRC (-7) variant (SEQ ID NO:144). In another embodiment, the variant is a MET variant of *P. tremuloides* isoprene synthase (SEQ ID NO:146). In another embodiment, the variant is a MET variant of *P. trichosarpa* isoprene synthase (SEQ ID NO:148). In another embodiment, the variant comprises one of more amino acid substitutions selected from the group consisting of V10M, F12S, T15A, E18G, V58I, V58F, L70Q, L70V, L70T, T71P, V79L, E89D, G94A, S119F, F120L, G127R, E175V, T212I, S257A, R262G, A266G, F280L, N297K, F305L, L319M, E323K, A328T, D342E, A359T, K366N, E368D, L374M, S396T, V418S, K438N, H440R, T442I, T442A, I449V, A469S, K500R, K505Q, G507S, S509N, F511Y, and N532K. In another embodiment, at least one amino acid substitution is a L70R substitution. In another embodiment, the variant comprises one of more amino acid substitutions selected from the group consisting of G127R/F511Y, L70Q/G94A/R262G/F305L, F12S/ T15A/E18G/N297K, S396T/T442I, V10M/E323K, F120L/A266G, K438N/K500R, V79L/S509N, E175V/S257A/E368D/A469S, T71P/L374M,

F280L/H440R, E89D/H440R, V58F/A328T/N532K, S119F/D342E/I449V, and K366N/G507S. In another embodiment, the polynucleotide sequence is contained within a plasmid. In another embodiment, the polynucleotide sequence is integrated into a chromosome of the host cell.

In another embodiment, the host is selected from the group consisting of gram-positive bacterial cells, gram-negative bacterial cells, filamentous fungal cells, and yeast cells. In another embodiment, the host is selected from the group consisting of *Escherichia sp.* (*E. coli*), *Pantoea sp.* (*P. citrea*), *Bacillus sp.* (*B. subtilis*), *Yarrowia sp.* (*Y. lipolytica*), and *Trichoderma* (*T. reesei*). In another embodiment, the host cell is cultured in a medium comprising a carbon source selected from the group consisting of glucose, glycerol, glycerine, dihydroxyacetone, yeast extract, biomass, molasses, sucrose, and oil. In another embodiment, the host cell further comprises a heterologous or native nucleic acid encoding an IDI polypeptide and/or a heterologous or native nucleic acid encoding a DXS polypeptide, optionally in combination with the native DXP pathway. In another embodiment, the host cell further comprises one or more nucleic acids encoding an IDI polypeptide and a DXS polypeptide. In another embodiment, the host cell comprises one vector encoding the isoprene synthase variant, the IDI polypeptide, and the DXS polypeptide. In another embodiment, the host cell further comprises a heterologous nucleic acid encoding an MVA pathway polypeptide selected from the group consisting of an MVA pathway polypeptide from *Saccharomyces cerevisia* and *Enterococcus faecalis*. In another embodiment, the host cell further comprises one or more nucleic acids encoding an MVA pathway polypeptide and a DXS polypeptide and wherein one vector encodes the isoprene synthase variant, the MVA pathway polypeptide, and the DXS polypeptide. In another embodiment, the host cell further comprises one or more nucleic acids encoding a DXS polypeptide, an IDI polypeptide, or one or more of the rest of the DXP pathway polypeptides, and a MVA pathway polypeptide.

[0019] In another aspect, the invention provides for methods of producing isoprene, comprising: (a) culturing the host cells comprising a heterologous polynucleotide sequence encoding an isoprene synthase variant in operable combination with a promoter, wherein the isoprene synthase variant comprises a substitution at a position corresponding to one or more residues (one, two, three, four, five, six, seven, eight, nine or ten) of a poplar isoprene synthase under suitable culture conditions for production of isoprene; and (b) producing the isoprene. In one embodiment, the method further comprises (c) recovering the isoprene. In another embodiment, the method further comprises (d) polymerizing isoprene.

[0020] In another aspect, the invention provides for methods of producing isoprene synthase, comprising: (a) providing: (i) a host cell; and (ii) a nucleic acid encoding an isoprene synthase variant in operable combination with a promoter, wherein the isoprene synthase variant comprises a substitution at a position corresponding to one or more residues (one, two, three, four, five, six, seven, eight, nine or ten) of a *P. alba* isoprene synthase of SEQ ID NO:120; (b) contacting the host cell with the nucleic acid to produce a transformed host cell; and (c) culturing the transformed host cells under suitable culture conditions for production of isoprene synthase.

BRIEF DESCRIPTION OF THE DRAWINGS

[0021] Figure 1 provides the coding sequence (SEQ ID NO:1) of kudzu (*Pueraria montana*) isoprene synthase, codon-optimized for expression in *Escherichia coli*.

[0022] Figure 2 provides the amino acid sequence (SEQ ID NO:2) of kudzu isoprene synthase.

[0023] Figure 3 provides the coding sequence (SEQ ID NO:6) of poplar (*Populus alba x tremula*) isoprene synthase, codon-optimized for expression in *Escherichia coli*.

[0024] Figure 4 provides the amino acid sequence (SEQ ID NO:7) of poplar (*Populus alba x tremula*) isoprene synthase.

[0025] Figure 5 provides a kudzu isoprene synthase homology model with the cysteine residues highlighted as space filling molecules.

[0026] Figure 6 provides a poplar isoprene synthase homology model with the cysteine residues highlighted as space filling molecules (dark grey). In addition, the cysteine residues from the kudzu model of Figure 5 are superimposed on the poplar model as space filling molecules (light grey).

[0027] Figure 7 provides a graph showing the growth curves of the kudzu IspS cysteine mutants of Example 5.

[0028] Figure 8 shows an SDS-PAGE analysis of kudzu IspS cysteine mutants from lysed cells. Pellet and supernatant fractions were prepared by centrifugation.

[0029] Figure 9 provides a map of plasmid MCM93 (pCR2.1-Kudzu).

[0030] Figure 10 provides the nucleotide sequence of plasmid MCM93 (SEQ ID NO:22).

[0031] Figure 11 provides a map of pET24D-Kudzu.

[0032] Figure 12 provides the nucleotide sequence of pET24D-Kudzu (SEQ ID NO:23).

[0033] Figure 13 shows an SDS-PAGE analysis of kudzu isoprene synthase-containing inclusion bodies. Lane M contains molecular weight markers, while the other lanes contain

increasing amounts of the purified inclusion body preparation. The kudzu isoprene synthase was estimated to have a purity of >90%.

[0034] Figure 14 provides graphs showing isoprene synthase activity of kudzu site evaluation library (SEL) members for positions Y22, A20 and S408. Most members show highly decreased activity, while conservative substitutions show a lesser decrease in activity. Panel A shows assay results for the Y22 library members in comparison with independent wild type samples (circled WT). Panel B shows assay results for the A20 library members in comparison to wild type samples (circled WT). Panel C shows assay results for the S408 library members, indicating that member S408D has 1.5 to 2-fold higher activity than the average of the wild type controls.

[0035] Figure 15 shows the MVA and DXP metabolic pathways for isoprene (based on F. Bouvier et al., Progress in Lipid Res. 44: 357-429, 2005). The following description includes alternative names for each polypeptide in the pathways and a reference that discloses an assay for measuring the activity of the indicated polypeptide (each of these references are each hereby incorporated by reference in their entireties, particularly with respect to assays for polypeptide activity for polypeptides in the MVA and DXP pathways). **Mevalonate Pathway:** **AACT**; Acetyl-CoA acetyltransferase, MvaE, EC 2.3.1.9. Assay: J. Bacteriol., 184: 2116–2122, 2002; **HMGS**; Hydroxymethylglutaryl-CoA synthase, MvaS, EC 2.3.3.10. Assay: J. Bacteriol., 184: 4065–4070, 2002; **HMGR**; 3-Hydroxy-3-methylglutaryl-CoA reductase, MvaE, EC 1.1.1.34. Assay: J. Bacteriol., 184: 2116–2122, 2002; **MVK**; Mevalonate kinase, ERG12, EC 2.7.1.36. Assay: Curr Genet 19:9-14, 1991. **PMK**; Phosphomevalonate kinase, ERG8, EC 2.7.4.2, Assay: Mol Cell Biol., 11:620–631, 1991; **DPMDC**; Diphosphomevalonate decarboxylase, MVD1, EC 4.1.1.33. Assay: Biochemistry, 33:13355-13362, 1994; **IDI**; Isopentenyl-diphosphate delta-isomerase, IDI1, EC 5.3.3.2. Assay: J. Biol. Chem. 264:19169-19175, 1989. **DXP Pathway:** **DXS**; 1-Deoxyxylulose-5-phosphate synthase, dxs, EC 2.2.1.7. Assay: PNAS, 94:12857-62, 1997; **DXR**; 1-Deoxy-D-xylulose 5-phosphate reductoisomerase, dxr, EC 2.2.1.7. Assay: Eur. J. Biochem. 269:4446–4457, 2002; **MCT**; 4-Diphosphocytidyl-2C-methyl-D-erythritol synthase, IspD, EC 2.7.7.60. Assay: PNAS, 97: 6451–6456, 2000; **CMK**; 4-Diphosphocytidyl-2-C-methyl-D-erythritol kinase, IspE, EC 2.7.1.148. Assay: PNAS, 97:1062-1067, 2000; **MCS**; 2C-Methyl-D-erythritol 2,4-cyclodiphosphate synthase, IspF, EC 4.6.1.12. Assay: PNAS, 96:11758-11763, 1999; **HDS**; 1-Hydroxy-2-methyl-2-(E)-butenyl 4-diphosphate synthase, ispG, EC 1.17.4.3. Assay: J. Org. Chem., 70:9168 -9174, 2005; **HDR**; 1-Hydroxy-2-methyl-2-(E)-butenyl 4-diphosphate reductase, IspH, EC 1.17.1.2. Assay: JACS, 126:12847-12855, 2004.

- [0036] Figure 16 provides a map of pDu27.
- [0037] Figure 17 provides the amino acid sequence (SEQ ID NO:118) of 6Xhis N-terminally tagged *P. alba* IspS in pDu27.
- [0038] Figure 18 provides the nucleotide sequence (SEQ ID NO:119) of plasmid pDu27.
- [0039] Figure 19 provides the amino acid sequence (SEQ ID NO:120) of full length *P. alba* IspS in pET24a. Underlined residues indicate the locations of N-terminal truncations in IspS in plasmids pDu39 through pDu43.
- [0040] Figure 20 provides the nucleotide sequence (SEQ ID NO:121) of plasmid *P. alba* pET24a.
- [0041] Figure 21 shows purified IspS displays a lower molecular weight “doublet” by SDS-PAGE analysis.
- [0042] Figure 22 shows tryptic peptides identified by mass spectrometry.
- [0043] Figure 23 provides maps of pDu40, pDu41, pDu42 and pDu43 harboring N-terminal truncations of *P. alba* IspS.
- [0044] Figure 24 provides a map of Pdu39 which is pET24a-*P. alba* MEA untagged (in strain MD09-173).
- [0045] Figure 25 provides the amino acid sequence (SEQ ID NO:122) of truncated “MEA” variant of *P. alba* IspS in pDu39.
- [0046] Figure 26 provides the nucleotide sequence (SEQ ID NO:123) of plasmid pDu39.
- [0047] Figure 27 provides the amino acid sequence (SEQ ID NO:124) of truncated “MSV” variant *P. alba* IspS in pDu41.
- [0048] Figure 28 provides the nucleotide sequence (SEQ ID NO:125) of plasmid pDu41 (pET24a-P.alba(MSV) Untagged).
- [0049] Figure 29 provides the amino acid sequence (SEQ ID NO:126) of truncated “MVS” variant *P. alba* IspS in pDu43.
- [0050] Figure 30 provides the nucleotide sequence (SEQ ID NO:127) of plasmid pDu43 (pET24a-P.alba(MVS) Untagged).
- [0051] Figure 31 provides the amino acid sequence (SEQ ID NO:128) of truncated “MTE” variant of *P. alba* IspS in pDu42.
- [0052] Figure 32 provides the nucleotide sequence (SEQ ID NO:129) of plasmid pDu42 (pET24a-P.alba(MTE) Untagged).
- [0053] Figure 33 provides the amino acid sequence (SEQ ID NO:130) of truncated “MNV” *P. alba* IspS in pDu40.

[0054] Figure 34 provides the nucleotide sequence (SEQ ID NO:131) of plasmid pDu40 (pET24a-P.alba(MNV) Untagged).

[0055] Figure 35 provides maps of MD09-161 and MD09-163, C-terminally TEV, 6XHis-tagged IspS variants.

[0056] Figure 36 provides the amino acid sequence (SEQ ID NO:132) of P. alba MEA(+)TEV in MD09-163.

[0057] Figure 37 provides the nucleotide sequence (SEQ ID NO:133) of plasmid MD09-163 (pET24a-P.alba MEA(+)TEV. CDS is underlined, TEV protease site is bold.

[0058] Figure 38 provides the amino acid sequence (SEQ ID NO:134) of P. alba FL(+)TEV in MD09-161.

[0059] Figure 39 provides the nucleotide sequence (SEQ ID NO:135) of plasmid MD09-161 (pET24a-P.alba FL(+)TEV. CDS is underlined, TEV protease site is bold.

[0060] Figure 40 shows a graph representing specific activities of MD09-167, Full length (FL), MD09-165, and Truncated isoprene synthase (MD09-173). Reactions were run at 30 °C for 15 minutes in a solution containing 100 mM Tris, 100 mM NaCl, 50 mM MgCl₂, 5 mM DMAPP, and 2.5-4.5 µg isoprene synthase in the supernatant of whole cell lysate.

[0061] Figure 41 shows graphs demonstrating Rate/[E] vs. [DMAPP]. X's represent MD09-173, circles represent MD09-167, diamonds represent MD09-165 and squares represent full length IspS.

[0062] Figure 42 shows a graph showing Isoprene synthase activity vs. [DMAPP]. X's represent data generated with MD09-173 truncated isoprene synthase. Circles represent data generated with MD09-167 isoprene synthase. Diamonds represent data generated with MD09-165 isoprene synthase. Squares represent data generated with full length isoprene synthase. Each data set was run in triplicate from independently grown cultures.

[0063] Figure 43 shows graphs demonstrating the effects of varying k_{cat} and K_M and K_i on reaction rate. In panel A, line 1 represents the rate equation of truncated isoprene synthase activity divided by the rate equation of full length isoprene synthase plotted at varying DMAPP concentrations. Line 2 represents the rate equation of full length isoprene synthase in which the k_{cat} has been substituted by the k_{cat} of the truncated isoprene synthase divided by the rate equation of the full length isoprene synthase. Line 3 represents the rate equation of full length isoprene synthase in which the K_M has been substituted by the K_M of the truncated isoprene synthase divided by the rate equation of the full length isoprene synthase. Line 4 represents the rate equation of full length isoprene synthase in which the K_i has been substituted by the K_i of the truncated isoprene synthase divided by the rate equation of the

full length isoprene synthase. Panel B shows a graph demonstrating data fit to the ratio of the rate equation of truncated isoprene synthase to full length isoprene synthase vs.[DMAPP].

[0064] Figure 44 shows a graph demonstrating growth inhibition of MCM531 by Mevalonic Acid (MVA). Cells were grown in TM3 medium in a microtiter plate with different concentrations of MVA. OD₆₀₀ of quadruplicate wells was measured at the indicated times.

[0065] Figure 45 shows graphs demonstrating DMAPP assays of L70 SSL plate. Dark bars represent either the full length (*P. alba* pET24a) or pDU39 (truncated) controls. The variants in wells C3 (27), D3 (39), or E3 (51) were chosen for further analysis.

[0066] Figure 46 shows a graph demonstrating the average specific activity of all variants selected for DMAPP assay with protein determination. Error bars show one standard deviation. All 3 L70R variants display higher activity than the control (WT).

[0067] Figure 47 shows a graph demonstrating the average specific activity of all 3 L70R variants compared to the “MEA” truncated *P. alba* IspS enzyme. Error bars show one standard deviation.

[0068] Figure 48 provides maps of plasmids pDu47-3, pDu47-4, and pDu47-5.

[0069] Figure 49 provides maps of plasmids pDu47-6, pDu47-7, and pDu48.

[0070] Figure 50 provides maps of plasmids pDu49, pDu50 and pDu50-4.

[0071] Figure 51 provides the amino acid sequence (SEQ ID NO:136) of *P. alba* TRC (-3) in pDu47-3.

[0072] Figure 52 provides the nucleotide sequence (SEQ ID NO:137) of plasmid pDu47-3.

[0073] Figure 53 provides the amino acid sequence (SEQ ID NO:138) of *P. alba* TRC (-4) in pDu47-4.

[0074] Figure 54 provides the nucleotide sequence (SEQ ID NO:139) of plasmid pDu47-4.

[0075] Figure 55 provides the amino acid sequence (SEQ ID NO:140) of *P. alba* TRC (-5) in pDu47-5.

[0076] Figure 56 provides the nucleotide sequence (SEQ ID NO:141) of plasmid pDu47-5.

[0077] Figure 57 provides the amino acid sequence (SEQ ID NO:142) of *P. alba* TRC (-6) in pDu47-6.

[0078] Figure 58 provides the nucleotide sequence (SEQ ID NO:143) of plasmid pDu47-6.

[0079] Figure 59 provides the amino acid sequence (SEQ ID NO:144) of *P. alba* TRC (-7) in pDu47-7.

[0080] Figure 60 provides the nucleotide sequence (SEQ ID NO:145) of plasmid pDu47-7.

- [0081] Figure 61 provides the amino acid sequence (SEQ ID NO:146) of *P. tremuloides* TRC (MET) in pDu48.
- [0082] Figure 62 provides the nucleotide sequence (SEQ ID NO:147) of plasmid pDu48.
- [0083] Figure 63 provides the amino acid sequence (SEQ ID NO:148) of *P. trichocarpa* (TRC) in pDu49.
- [0084] Figure 64 provides the nucleotide sequence (SEQ ID NO:149) of plasmid pDu49.
- [0085] Figure 65 provides the amino acid sequence (SEQ ID NO:150) of Kudzu TRC (MEA) in pDu50.
- [0086] Figure 66 provides the nucleotide sequence (SEQ ID NO:151) of plasmid pDu50.
- [0087] Figure 67 provides the amino acid sequence (SEQ ID NO:152) of KudzuTRC (-4) in pDu50-4.
- [0088] Figure 68 provides the nucleotide sequence (SEQ ID NO:153) of plasmid pDu50-4.
- [0089] Figure 69 shows graphs demonstrating raw and OD-normalized data from DMAPP assay of truncated variants of IspS.
- [0090] Figure 70 shows a graph representing the specific activity of IspS truncations. *P. alba*, *P. tremuloides* and *P. trichocarpa* truncations were compared for specific activity relative to the *P. alba* "full length" variant.
- [0091] Figure 71 provides a map of plasmid p9795.
- [0092] Figure 72 provides the nucleotide sequence of plasmid p9795 (SEQ ID NO:154).
- [0093] Figure 73 provides a map of plasmid pTrcKudzu.
- [0094] Figure 74 provides the nucleotide sequence (SEQ ID NO:155) of plasmid pTrcKudzu.
- [0095] Figure 75 provides a map of plasmid pMAL-C4X.
- [0096] Figure 76 provides the nucleotide sequence (SEQ ID NO:156) of plasmid pMAL-C4X.
- [0097] Figure 77 provides a map of plasmid pMAL-C4X-Kudzu.
- [0098] Figure 78 provides the nucleotide sequence (SEQ ID NO:157) of plasmid pMAL-C4X-Kudzu.
- [0099] Figure 79 provides maps of plasmids pET24 *P. tremuloides* pET24a and *P. trichocarpa* pET24a.
- [0100] Figure 80 provides the amino acid sequence (SEQ ID NO:158) of *P. tremuloides* IspS in *P. trichocarpa* pET24a.
- [0101] Figure 81 provides the nucleotide sequence (SEQ ID NO:159) of plasmid *P. tremuloides* pET24a.

[0102] Figure 82 provides the amino acid sequence (SEQ ID NO:160) of *P. trichocharpa* IspS in *P. trichocharpa* pET24a.

[0103] Figure 83 provides the nucleotide sequence (SEQ ID NO:161) of plasmid *P. trichocharpa* pET24a.

[0104] Figure 84 provides maps of plasmids pDu30, pDu31, and pDu32.

[0105] Figure 85 provides the amino acid sequence (SEQ ID NO:162) of IspS variant *P.alba*TRC-pET200 in pDu30.

[0106] Figure 86 provides the nucleotide sequence (SEQ ID NO:163) of pDu30.

[0107] Figure 87 provides the amino acid sequence (SEQ ID NO:164) of IspS variant *P.trem*TRC-pET200 in pDu31.

[0108] Figure 88 provides the nucleotide sequence (SEQ ID NO:165) of pDu31.

[00100] Figure 89 provides the amino acid sequence (SEQ ID NO:166) of IspS variant *P.trich*TRC-pET200 in pDu32.

[0109] Figure 90 provides the nucleotide sequence (SEQ ID NO:167) of pDu32.

[0110] Figure 91 provides the three-dimensional structure of *P. tremuloides* IspS shown as a dimer. Chain A is in dark gray, chain B is in medium gray and the single magnesium ion found in each active site is light gray.

[0111] Figure 92 provides a monomer view of the structure of *P. tremuloides* IspS. The magnesium is shown as a light gray sphere and the N- and C-terminals are indicated.

[0112] Figure 93 shows the structural alignments between (A) BdpS and LS, (B) BdpS and poplar IspS, (C) LS and poplar IspS, and (D) TEAS and poplar IspS. In each case the first structure is in light gray and the second is in dark gray. Divalent cations are shown as spheres.

[0113] Figure 94 shows the three dimensional structure of loops in BdpS and LS. Panel A shows the N-terminal loop of Ls in light gray and the N-terminal loop of BdpS in dark gray. Panel B shows that Loop I and Loop II are structurally homologous.

[0114] Figure 95 shows the N-terminal loop of BdpS (dark gray) and poplar IspS (light gray) are structurally divergent. Panel A shows the N-terminal loop and panel B shows Loop I and Loop II.

[0115] Figure 96 shows the N-terminal loop of LS (light gray) and poplar IspS (dark gray) are structurally divergent. Panel A shows the N-terminal loop and panel B shows Loop I and Loop II.

[0116] Figure 97 shows the N-terminal loop of TEAS (light gray) and poplar IspS (dark gray) are structurally divergent. Panel A shows the N-terminal loop and panel B shows Loop I and Loop II. Loop I is disordered in TEAS.

GENERAL DESCRIPTION OF THE INVENTION

[0117] The present invention provides methods and compositions comprising at least one isoprene synthase enzyme with improved catalytic activity and/or solubility. In particular, the present invention provides variant plant isoprene synthases for increased isoprene production in microbial host cells. Biosynthetically produced isoprene of the present invention finds use in the manufacture of rubber and elastomers.

[0118] Unless otherwise indicated, the practice of the present invention involves conventional techniques commonly used in molecular biology, microbiology, and recombinant DNA, which are within the skill of the art. Such techniques are known to those of skill in the art and are described in numerous texts and reference works (See e.g., Sambrook et al., "Molecular Cloning: A Laboratory Manual," Second Edition, Cold Spring Harbor, 1989; and Ausubel et al., "Current Protocols in Molecular Biology," 1987).

[0119] Unless defined otherwise herein, all technical and scientific terms used herein have the same meaning as commonly understood by one of ordinary skill in the art to which this invention pertains. For example, Singleton and Sainsbury, Dictionary of Microbiology and Molecular Biology, 2d Ed., John Wiley and Sons, NY (1994); and Hale and Marham, The Harper Collins Dictionary of Biology, Harper Perennial, NY (1991) provide those of skill in the art with a general dictionaries of many of the terms used in the invention. Although any methods and materials similar or equivalent to those described herein find use in the practice of the present invention, the preferred methods and materials are described herein.

Accordingly, the terms defined immediately below are more fully described by reference to the Specification as a whole.

[0120] Furthermore, the headings provided herein are not limitations of the various aspects or embodiments of the invention, which can be had by reference to the specification as a whole. Accordingly, the terms defined immediately below are more fully defined by reference to the specification as a whole. Nonetheless, in order to facilitate understanding of the invention, a number of terms are defined below.

Definitions

[0121] Unless defined otherwise herein, all technical and scientific terms used herein have the same meaning as commonly understood by one of ordinary skill in the art to which this invention pertains. Although any methods and materials similar or equivalent to those described herein find use in the practice of the present invention, the preferred methods and materials are described herein. Accordingly, the terms defined immediately below are more fully described by reference to the Specification as a whole.

[0122] As used herein, the singular terms “a,” “an,” and “the” include the plural reference unless the context clearly indicates otherwise. Unless otherwise indicated, nucleic acids are written left to right in 5' to 3' orientation; amino acid sequences are written left to right in amino to carboxy orientation, respectively. It is to be understood that this invention is not limited to the particular methodology, protocols, and reagents described, as these may vary, depending upon the context they are used by those of skill in the art.

[0123] It is intended that every maximum numerical limitation given throughout this specification includes every lower numerical limitation, as if such lower numerical limitations were expressly written herein. Every minimum numerical limitation given throughout this specification will include every higher numerical limitation, as if such higher numerical limitations were expressly written herein. Every numerical range given throughout this specification will include every narrower numerical range that falls within such broader numerical range, as if such narrower numerical ranges were all expressly written herein.

[0124] All documents cited are, in relevant part, incorporated herein by reference. However, the citation of any document is not to be construed as an admission that it is prior art with respect to the present invention.

[0125] As used herein, the term 2-methyl-1,3-butadiene (CAS# 78-79-5) (“isoprene”) refers to the direct and final volatile C5 hydrocarbon product from the elimination of pyrophosphate from 3,3-dimethylallyl pyrophosphate (DMAPP), and does not involve the linking or polymerization of [an] IPP molecule(s) to [a] DMAPP molecule(s). As used herein, the terms “isoprene synthase,” and “IspS,” refer to the enzymes that catalyze the elimination of pyrophosphate from diemethylallyl diphosphate (DMAPP) to form isoprene. In some preferred embodiments, the IspS is an enzyme obtained from plants such as kudzu, poplar or red oak. In some embodiments, the term “IspS” refers to a naturally occurring mature enzyme or portion thereof.

[0126] Related (and derivative) proteins comprise “variant proteins.” In some preferred embodiments, variant proteins differ from a parent protein (e.g., kudzu IspS set forth as SEQ ID NO:2 or poplar IspS) and one another by a small number of amino acid residues. The number of differing amino acid residues may be one or more, preferably 1, 2, 3, 4, 5, 10, 15, 20, 30, 40, 50, or more amino acid residues. In some preferred embodiments, the number of different amino acids between variants is between 1 and 10. In some particularly preferred embodiments, related proteins and particularly variant proteins comprise at least 35%, 40%, 45%, 50%, 55%, 60%, 65%, 70%, 75%, 80%, 85%, 90%, 95%, 97%, 98%, or 99% amino acid sequence identity. Additionally, a related protein or a variant protein as used herein

refers to a protein that differs from another related protein or a parent protein in the number of prominent regions. For example, in some embodiments, variant proteins have 1, 2, 3, 4, 5, or 10 corresponding prominent regions that differ from the parent protein.

[0127] Several methods are known in the art that are suitable for generating variants of the enzymes of the present invention, including but not limited to site-saturation mutagenesis, scanning mutagenesis, insertional mutagenesis, random mutagenesis, site-directed mutagenesis, and directed-evolution, as well as various other recombinatorial approaches.

[0128] Characterization of wild-type and mutant proteins is accomplished via any means or “test” suitable and is preferably based on the assessment of properties of interest. For example one or more of the following properties are assessed in some embodiments of the present invention: pH stability; temperature stability; oxidative stability; proteolytic stability; solubility; K_m and/or specific activity of the conversion of DMAPP to isoprene *in vitro*; K_m and/or specific activity of the conversion of DMAPP to isoprene *in vivo* in the context of a host organism (e.g., *E. coli*); and expression of enzyme(s) of the DXP pathway and/or the MVA pathway. Indeed, it is contemplated that enzymes having various degrees of stability, solubility, activity, and/or expression level in one or more of test conditions will find use in the present invention.

[0129] As used herein the term “gene” refers to a polynucleotide (e.g., a DNA segment) that encodes a polypeptide and includes regions preceding and following the coding regions as well as intervening sequences (introns) between individual coding segments (exons).

[0130] As used herein, “homologous genes” refers to a pair of genes from different, but usually related species, which correspond to each other and which are identical or very similar to each other. The term encompasses genes that are separated by speciation (*i.e.*, the development of new species) (e.g., orthologous genes), as well as genes that have been separated by genetic duplication (e.g., paralogous genes).

[0131] As used herein, “ortholog” and “orthologous genes” refer to genes in different species that have evolved from a common ancestral gene (*i.e.*, a homologous gene) by speciation. Typically, orthologs retain the same function during the course of evolution. Identification of orthologs finds use in the reliable prediction of gene function in newly sequenced genomes.

[0132] As used herein, “paralog” and “paralogous genes” refer to genes that are related by duplication within a genome. While orthologs retain the same function through the course of evolution, paralogs evolve new functions, even though some functions are often related to the original one.

[0133] As used herein, “homology” refers to sequence similarity or identity, with identity being preferred. This homology is determined using standard techniques known in the art (*See e.g.*, Smith and Waterman, *Adv Appl Math*, 2:482, 1981; Needleman and Wunsch, *J Mol Biol*, 48:443, 1970; Pearson and Lipman, *Proc Natl Acad Sci USA*, 85:2444, 1988; programs such as GAP, BESTFIT, FASTA, and TFASTA in the Wisconsin Genetics Software Package, Genetics Computer Group, Madison, WI; and Devereux *et al.*, *Nucl Acid Res*, 12:387-395, 1984).

[0134] As used herein, an “analogous sequence” is one wherein the function of the gene is essentially the same as the gene based on the kudzu isoprene synthase (IspS) or poplar IspS (IspS). Additionally, analogous genes include at least 45%, 50%, 55%, 60%, 65%, 70%, 75%, 80%, 85%, 90%, 95%, 97%, 98%, 99% or 100% sequence identity with the sequence of the kudzu isoprene synthase. In additional embodiments more than one of the above properties applies to the sequence. Analogous sequences are determined by known methods of sequence alignment. A commonly used alignment method is BLAST, although as indicated above and below, there are other methods that also find use in aligning sequences.

[0135] One example of a useful algorithm is PILEUP. PILEUP creates a multiple sequence alignment from a group of related sequences using progressive, pair-wise alignments. It can also plot a tree showing the clustering relationships used to create the alignment. PILEUP uses a simplification of the progressive alignment method of Feng and Doolittle (Feng and Doolittle, *J Mol Evol*, 35:351-360, 1987). The method is similar to that described by Higgins and Sharp (Higgins and Sharp, *CABIOS* 5:151-153, 1989). Useful PILEUP parameters including a default gap weight of 3.00, a default gap length weight of 0.10, and weighted end gaps.

[0136] Another example of a useful algorithm is the BLAST algorithm, described by Altschul *et al.*, (Altschul *et al.*, *J Mol Biol*, 215:403-410, 1990; and Karlin *et al.*, *Proc Natl Acad Sci USA*, 90:5873-5787, 1993). A particularly useful BLAST program is the WU-BLAST-2 program (*See*, Altschul *et al.*, *Meth Enzymol*, 266:460-480, 1996). WU-BLAST-2 uses several search parameters, most of which are set to the default values. The adjustable parameters are set with the following values: overlap span = 1, overlap fraction = 0.125, word threshold (T) = 11. The HSP S and HSP S2 parameters are dynamic values and are established by the program itself depending upon the composition of the particular sequence and composition of the particular database against which the sequence of interest is being searched. However, the values may be adjusted to increase sensitivity. A % amino acid sequence identity value is determined by the number of matching identical residues divided

by the total number of residues of the “longer” sequence in the aligned region. The “longer” sequence is the one having the most actual residues in the aligned region (gaps introduced by WU-Blast-2 to maximize the alignment score are ignored).

[0137] Thus, “percent (%) nucleic acid sequence identity” is defined as the percentage of nucleotide residues in a candidate sequence that are identical to the nucleotide residues of the starting sequence (*i.e.*, the sequence of interest). A preferred method utilizes the BLASTN module of WU-BLAST-2 set to the default parameters, with overlap span and overlap fraction set to 1 and 0.125, respectively.

[0138] As used herein, the term “hybridization” refers to the process by which a strand of nucleic acid joins with a complementary strand through base pairing, as known in the art.

[0139] A nucleic acid sequence is considered to be “selectively hybridizable” to a reference nucleic acid sequence if the two sequences specifically hybridize to one another under moderate to high stringency hybridization and wash conditions. Hybridization conditions are based on the melting temperature (T_m) of the nucleic acid binding complex or probe. For example, “maximum stringency” typically occurs at about $T_m-5^{\circ}\text{C}$ (5° below the T_m of the probe); “high stringency” at about $5-10^{\circ}\text{C}$ below the T_m ; “intermediate stringency” at about $10-20^{\circ}\text{C}$ below the T_m of the probe; and “low stringency” at about $20-25^{\circ}\text{C}$ below the T_m . Functionally, maximum stringency conditions may be used to identify sequences having strict identity or near-strict identity with the hybridization probe; while intermediate or low stringency hybridization can be used to identify or detect polynucleotide sequence homologs.

[0140] Moderate and high stringency hybridization conditions are well known in the art. An example of high stringency conditions includes hybridization at about 42°C in 50% formamide, 5X SSC, 5X Denhardt’s solution, 0.5% SDS and 100 $\mu\text{g/ml}$ denatured carrier DNA followed by washing two times in 2X SSC and 0.5% SDS at room temperature and two additional times in 0.1X SSC and 0.5% SDS at 42°C . An example of moderate stringent conditions include an overnight incubation at 37°C in a solution comprising 20% formamide, 5 x SSC (150mM NaCl, 15 mM trisodium citrate), 50 mM sodium phosphate (pH 7.6), 5 x Denhardt’s solution, 10% dextran sulfate and 20 mg/ml denatured sheared salmon sperm DNA, followed by washing the filters in 1x SSC at about $37 - 50^{\circ}\text{C}$. Those of skill in the art know how to adjust the temperature, ionic strength, etc. as necessary to accommodate factors such as probe length and the like.

[0141] As used herein, “recombinant” includes reference to a cell or vector, that has been modified by the introduction of a heterologous nucleic acid sequence or that the cell is

derived from a cell so modified. Thus, for example, recombinant cells express genes that are not found in identical form within the native (non-recombinant) form of the cell or express native genes that are otherwise abnormally expressed, under expressed or not expressed at all as a result of deliberate human intervention. "Recombination," "recombining," and generating a "recombined" nucleic acid are generally the assembly of two or more nucleic acid fragments wherein the assembly gives rise to a chimeric gene.

[0142] In a preferred embodiment, mutant DNA sequences are generated with site saturation mutagenesis in at least one codon. In another preferred embodiment, site saturation mutagenesis is performed for two or more codons. In a further embodiment, mutant DNA sequences have more than 50%, more than 55%, more than 60%, more than 65%, more than 70%, more than 75%, more than 80%, more than 85%, more than 90%, more than 95%, or more than 98% homology with the wild-type sequence. In alternative embodiments, mutant DNA is generated *in vivo* using any known mutagenic procedure such as, for example, radiation, nitrosoguanidine and the like. The desired DNA sequence is then isolated and used in the methods provided herein.

[0143] As used herein, the term "target sequence" refers to a DNA sequence in the host cell that encodes the sequence where it is desired for the incoming sequence to be inserted into the host cell genome. In some embodiments, the target sequence encodes a functional wild-type gene or operon, while in other embodiments the target sequence encodes a functional mutant gene or operon, or a non-functional gene or operon.

[0144] As used herein, a "flanking sequence" refers to any sequence that is either upstream or downstream of the sequence being discussed (*e.g.*, for genes A-B-C, gene B is flanked by the A and C gene sequences). In a preferred embodiment, the incoming sequence is flanked by a homology box on each side. In another embodiment, the incoming sequence and the homology boxes comprise a unit that is flanked by stuffer sequence on each side. In some embodiments, a flanking sequence is present on only a single side (either 3' or 5'), but in preferred embodiments, it is on each side of the sequence being flanked. In some embodiments, a flanking sequence is present on only a single side (either 3' or 5'), while in preferred embodiments, it is present on each side of the sequence being flanked.

[0145] As used herein, the term "stuffer sequence" refers to any extra DNA that flanks homology boxes (typically vector sequences). However, the term encompasses any non-homologous DNA sequence. Not to be limited by any theory, a stuffer sequence provides a noncritical target for a cell to initiate DNA uptake.

[0146] As used herein, the terms “amplification” and “gene amplification” refer to a process by which specific DNA sequences are disproportionately replicated such that the amplified gene becomes present in a higher copy number than was initially present in the genome. In some embodiments, selection of cells by growth in the presence of a drug (*e.g.*, an inhibitor of an inhibitable enzyme) results in the amplification of either the endogenous gene encoding the gene product required for growth in the presence of the drug or by amplification of exogenous (*i.e.*, input) sequences encoding this gene product, or both.

[0147] “Amplification” is a special case of nucleic acid replication involving template specificity. It is to be contrasted with non-specific template replication (*i.e.*, replication that is template-dependent but not dependent on a specific template). Template specificity is here distinguished from fidelity of replication (*i.e.*, synthesis of the proper polynucleotide sequence) and nucleotide (ribo- or deoxyribo-) specificity. Template specificity is frequently described in terms of “target” specificity. Target sequences are “targets” in the sense that they are sought to be sorted out from other nucleic acid. Amplification techniques have been designed primarily for this sorting out.

[0148] As used herein, the term “co-amplification” refers to the introduction into a single cell of an amplifiable marker in conjunction with other gene sequences (*i.e.*, comprising one or more non-selectable genes such as those contained within an expression vector) and the application of appropriate selective pressure such that the cell amplifies both the amplifiable marker and the other, non-selectable gene sequences. The amplifiable marker may be physically linked to the other gene sequences or alternatively two separate pieces of DNA, one containing the amplifiable marker and the other containing the non-selectable marker, may be introduced into the same cell.

[0149] As used herein, the terms “amplifiable marker,” “amplifiable gene,” and “amplification vector” refer to a gene or a vector encoding a gene, which permits the amplification of that gene under appropriate growth conditions.

[0150] “Template specificity” is achieved in most amplification techniques by the choice of enzyme. Amplification enzymes are enzymes that, under conditions they are used, will process only specific sequences of nucleic acid in a heterogeneous mixture of nucleic acid. For example, in the case of Q β replicase, MDV-1 RNA is the specific template for the replicase (*See e.g.*, Kacian *et al.*, Proc Natl Acad Sci USA 69:3038, 1972) and other nucleic acids are not replicated by this amplification enzyme. Similarly, in the case of T7 RNA polymerase, this amplification enzyme has a stringent specificity for its own promoters (*See*, Chamberlin *et al.*, Nature 228:227, 1970). In the case of T4 DNA ligase, the enzyme will not

ligate the two oligonucleotides or polynucleotides, where there is a mismatch between the oligonucleotide or polynucleotide substrate and the template at the ligation junction (*See*, Wu and Wallace, *Genomics* 4:560, 1989). Finally, *Taq* and *Pfu* polymerases, by virtue of their ability to function at high temperature, are found to display high specificity for the sequences bounded and thus defined by the primers; the high temperature results in thermodynamic conditions that favor primer hybridization with the target sequences and not hybridization with non-target sequences.

[0151] As used herein, the term “amplifiable nucleic acid” refers to nucleic acids, which may be amplified by any amplification method. It is contemplated that “amplifiable nucleic acid” will usually comprise “sample template.”

[0152] As used herein, the term “sample template” refers to nucleic acid originating from a sample, which is analyzed for the presence of “target” (defined below). In contrast, “background template” is used in reference to nucleic acid other than sample template, which may or may not be present in a sample. Background template is most often inadvertent. It may be the result of carryover, or it may be due to the presence of nucleic acid contaminants sought to be purified away from the sample. For example, nucleic acids from organisms other than those to be detected may be present as background in a test sample.

[0153] As used herein, the term “primer” refers to an oligonucleotide, whether occurring naturally as in a purified restriction digest or produced synthetically, which is capable of acting as a point of initiation of synthesis when placed under conditions in which synthesis of a primer extension product which is complementary to a nucleic acid strand is induced, (*i.e.*, in the presence of nucleotides and an inducing agent such as DNA polymerase and at a suitable temperature and pH). The primer is preferably single stranded for maximum efficiency in amplification, but may alternatively be double stranded. If double stranded, the primer is first treated to separate its strands before being used to prepare extension products. Preferably, the primer is an oligodeoxyribonucleotide. The primer must be sufficiently long to prime the synthesis of extension products in the presence of the inducing agent. The exact lengths of the primers will depend on many factors, including temperature, source of primer and the use of the method.

[0154] As used herein, the term “probe” refers to an oligonucleotide (*i.e.*, a sequence of nucleotides), whether occurring naturally as in a purified restriction digest or produced synthetically, recombinantly or by PCR amplification, which is capable of hybridizing to another oligonucleotide of interest. A probe may be single-stranded or double-stranded. Probes are useful in the detection, identification and isolation of particular gene sequences. It

is contemplated that any probe used in the present invention will be labeled with any “reporter molecule,” so that is detectable in any detection system, including, but not limited to enzyme (*e.g.*, ELISA, as well as enzyme-based histochemical assays), fluorescent, radioactive, and luminescent systems. It is not intended that the present invention be limited to any particular detection system or label.

[0155] As used herein, the term “target,” when used in reference to the polymerase chain reaction, refers to the region of nucleic acid bounded by the primers used for polymerase chain reaction. Thus, the “target” is sought to be sorted out from other nucleic acid sequences. A “segment” is defined as a region of nucleic acid within the target sequence.

[0156] As used herein, in one embodiment, the term “polymerase chain reaction” (“PCR”) refers to the methods of U.S. Patent Nos. 4,683,195 4,683,202, and 4,965,188, hereby incorporated by reference, which include methods for increasing the concentration of a segment of a target sequence in a mixture of genomic DNA without cloning or purification. This process for amplifying the target sequence consists of introducing a large excess of two oligonucleotide primers to the DNA mixture containing the desired target sequence, followed by a precise sequence of thermal cycling in the presence of a DNA polymerase. The two primers are complementary to their respective strands of the double stranded target sequence. To effect amplification, the mixture is denatured and the primers then annealed to their complementary sequences within the target molecule. Following annealing, the primers are extended with a polymerase so as to form a new pair of complementary strands. The steps of denaturation, primer annealing and polymerase extension can be repeated many times (*i.e.*, denaturation, annealing and extension constitute one “cycle”; there can be numerous “cycles”) to obtain a high concentration of an amplified segment of the desired target sequence. The length of the amplified segment of the desired target sequence is determined by the relative positions of the primers with respect to each other, and therefore, this length is a controllable parameter. By virtue of the repeating aspect of the process, the method is referred to as the “polymerase chain reaction” (hereinafter “PCR”). Because the desired amplified segments of the target sequence become the predominant sequences (in terms of concentration) in the mixture, they are said to be “PCR amplified”.

[0157] As used herein, the term “amplification reagents” refers to those reagents (deoxyribonucleotide triphosphates, buffer, etc.), needed for amplification except for primers, nucleic acid template and the amplification enzyme. Typically, amplification reagents along with other reaction components are placed and contained in a reaction vessel (test tube, microwell, etc.).

[0158] With PCR, it is possible to amplify a single copy of a specific target sequence in genomic DNA to a level detectable by several different methodologies (*e.g.*, hybridization with a labeled probe; incorporation of biotinylated primers followed by avidin-enzyme conjugate detection; incorporation of ³²P-labeled deoxynucleotide triphosphates, such as dCTP or dATP, into the amplified segment). In addition to genomic DNA, any oligonucleotide or polynucleotide sequence can be amplified with the appropriate set of primer molecules. In particular, the amplified segments created by the PCR process itself are, themselves, efficient templates for subsequent PCR amplifications.

[0159] As used herein, the terms “PCR product,” “PCR fragment,” and “amplification product” refer to the resultant mixture of compounds after two or more cycles of the PCR steps of denaturation, annealing and extension are complete. These terms encompass the case where there has been amplification of one or more segments of one or more target sequences.

[0160] As used herein, the term “RT-PCR” refers to the replication and amplification of RNA sequences. In this method, reverse transcription is coupled to PCR, most often using a one enzyme procedure in which a thermostable polymerase is employed, as described in U.S. Patent No. 5,322,770, herein incorporated by reference. In RT-PCR, the RNA template is converted to cDNA due to the reverse transcriptase activity of the polymerase, and then amplified using the polymerizing activity of the polymerase (*i.e.*, as in other PCR methods).

[0161] As used herein, the terms “restriction endonucleases” and “restriction enzymes” refer to bacterial enzymes, each of which cut double-stranded DNA at or near a specific nucleotide sequence.

[0162] A “restriction site” refers to a nucleotide sequence recognized and cleaved by a given restriction endonuclease and is frequently the site for insertion of DNA fragments. In certain embodiments of the invention restriction sites are engineered into the selective marker and into 5' and 3' ends of the DNA construct.

[0163] As used herein, the term “chromosomal integration” refers to the process whereby an incoming sequence is introduced into the chromosome of a host cell. The homologous regions of the transforming DNA align with homologous regions of the chromosome. Subsequently, the sequence between the homology boxes is replaced by the incoming sequence in a double crossover (*i.e.*, homologous recombination). In some embodiments of the present invention, homologous sections of an inactivating chromosomal segment of a DNA construct align with the flanking homologous regions of the indigenous chromosomal region of the *Escherichia* chromosome. Subsequently, the indigenous chromosomal region is deleted by the DNA construct in a double crossover (*i.e.*, homologous recombination).

[0164] “Homologous recombination” means the exchange of DNA fragments between two DNA molecules or paired chromosomes at the site of identical or nearly identical nucleotide sequences. In a preferred embodiment, chromosomal integration is homologous recombination. “Homologous sequences” as used herein means a nucleic acid or polypeptide sequence having 100%, 99%, 98%, 97%, 96%, 95%, 94%, 93%, 92%, 91%, 90%, 88%, 85%, 80%, 75%, or 70% sequence identity to another nucleic acid or polypeptide sequence when optimally aligned for comparison. In some embodiments, homologous sequences have between 85% and 100% sequence identity, while in other embodiments there is between 90% and 100% sequence identity, and in more preferred embodiments, there is 95% and 100% sequence identity.

[0165] As used herein “amino acid” refers to peptide or protein sequences or portions thereof. The terms “protein,” “peptide,” and “polypeptide” are used interchangeably.

[0166] As used herein, the term “heterologous protein” refers to a protein or polypeptide that does not naturally occur in the host cell. Examples of heterologous proteins include enzymes such as isoprene synthases. In some embodiments, the genes encoding the proteins are naturally occurring genes, while in other embodiments mutated and/or synthetic genes are used.

[0167] As used herein, “homologous protein” refers to a protein or polypeptide native or naturally occurring in a cell. In preferred embodiments, the cell is a Gram-negative cell, while in particularly preferred embodiments the cell is an *Escherichia* host cell.

[0168] An enzyme is “overexpressed” in a host cell if the enzyme is expressed in the cell at a higher level than the level at which it is expressed in a corresponding wild-type cell.

[0169] The terms “protein” and “polypeptide” are used interchangeably herein. The 3-letter code for amino acids as defined in conformity with the IUPAC-IUB Joint Commission on Biochemical Nomenclature (JCBN) is used throughout this disclosure. It is also understood that a polypeptide may be coded for by more than one nucleotide sequence due to the degeneracy of the genetic code.

[0170] The term “mature” form of a protein or peptide refers to the final functional form of the protein or peptide. To exemplify, a mature form of kudzu isoprene synthase includes the amino acid sequence of SEQ ID NO:2.

[0171] The term “precursor” form of a protein or peptide refers to a mature form of the protein having a prosequence operably linked to the amino or carboxyl terminus of the protein. The precursor may also have a “signal” sequence operably linked, to the amino terminus of the prosequence. The precursor may also have additional polynucleotides that are

involved in post-translational activity (*e.g.*, polynucleotides cleaved therefrom to leave the mature form of a protein or peptide).

[0172] “Naturally occurring enzyme” refers to an enzyme having the unmodified amino acid sequence identical to that found in nature. Naturally occurring enzymes include native enzymes, those enzymes naturally expressed or found in the particular microorganism.

[0173] The term “identical” in the context of two nucleic acids or polypeptide sequences refers to the residues in the two sequences that are the same when aligned for maximum correspondence, as measured using one of the following sequence comparison or analysis algorithms.

[0174] The term “optimal alignment” refers to the alignment giving the highest percent identity score.

[0175] “Percent sequence identity,” “percent amino acid sequence identity,” “percent gene sequence identity,” and/or “percent nucleic acid/polynucleotide sequence identity,” with respect to two amino acid, polynucleotide and/or gene sequences (as appropriate), refer to the percentage of residues that are identical in the two sequences when the sequences are optimally aligned. Thus, 80% amino acid sequence identity means that 80% of the amino acids in two optimally aligned polypeptide sequences are identical.

[0176] The phrase “substantially identical” in the context of two nucleic acids or polypeptides thus refers to a polynucleotide or polypeptide that comprising at least 70% sequence identity, preferably at least 75%, preferably at least 80%, preferably at least 85%, preferably at least 90%, preferably at least 95%, preferably at least 97% , preferably at least 98% and preferably at least 99% sequence identity as compared to a reference sequence using the programs or algorithms (*e.g.*, BLAST, ALIGN, CLUSTAL) using standard parameters. One indication that two polypeptides are substantially identical is that the first polypeptide is immunologically cross-reactive with the second polypeptide. Typically, polypeptides that differ by conservative amino acid substitutions are immunologically cross-reactive. Thus, a polypeptide is substantially identical to a second polypeptide, for example, where the two peptides differ only by a conservative substitution. Another indication that two nucleic acid sequences are substantially identical is that the two molecules hybridize to each other under stringent conditions (*e.g.*, within a range of medium to high stringency).

[0177] The term “isolated” or “purified” refers to a material that is removed from its original environment (*e.g.*, the natural environment if it is naturally occurring). For example, the material is said to be “purified” when it is present in a particular composition in a higher or lower concentration than exists in a naturally occurring or wild type organism (*e.g.*, kudzu)

or in combination with components not normally present upon expression from a naturally occurring or wild type organism. For example, a naturally-occurring polynucleotide or polypeptide present in a living animal is not isolated, but the same polynucleotide or polypeptide, separated from some or all of the coexisting materials in the natural system, is isolated. Such polynucleotides could be part of a vector, and/or such polynucleotides or polypeptides could be part of a composition, and still be isolated in that such vector or composition is not part of its natural environment. In preferred embodiments, a nucleic acid or protein is said to be purified, for example, if it gives rise to essentially one band in an electrophoretic gel or blot.

[0178] The term “isolated”, when used in reference to a DNA sequence, refers to a DNA sequence that has been removed from its natural genetic milieu and is thus free of other extraneous or unwanted coding sequences, and is in a form suitable for use within genetically engineered protein production systems. Similarly, the term “isolated”, when used in reference to a recombinant DNA sequence, refers to a DNA sequence that has been removed from the genetic milieu of the host organism and is thus free of other extraneous or unwanted coding sequences (e.g., kudzu IspS expression vector propagated in *E. coli*). Such isolated molecules are those that are separated from their natural environment and include cDNA and genomic clones. Isolated DNA molecules of the present invention are free of other genes with which they are ordinarily associated, but may include naturally occurring 5' and 3' untranslated regions such as promoters and terminators. The identification of associated regions will be evident to one of ordinary skill in the art (*See e.g.*, Dynan and Tijan, *Nature* 316:774-78, 1985). The term “an isolated DNA sequence” is alternatively referred to as “a cloned DNA sequence”.

[0179] The term “isolated,” when used in reference to a protein, refers to a protein that is found in a condition other than its native environment. In a preferred form, the isolated protein is substantially free of other proteins, particularly other homologous proteins. Similarly, the term “isolated”, when used in reference to a recombinantly produced protein, refers to a protein that has been removed from the proteinaceous milieu of the host organism and is thus free of other extraneous or unwanted proteins (e.g., recombinant kudzu IspS produced in *E. coli*). An isolated protein is more than 10% pure, preferably more than 20% pure, and even more preferably more than 30% pure, as determined by SDS-PAGE. Further aspects of the invention encompass the protein in a highly purified form (*i.e.*, more than 40% pure, more than 60% pure, more than 80% pure, more than 90% pure, more than 95% pure, more than 97% pure, and even more than 99% pure), as determined by SDS-PAGE.

[0180] The following cassette mutagenesis method may be used to facilitate the construction of the enzyme variants of the present invention, although other methods may be used. First, as described herein, a naturally-occurring gene encoding the enzyme is obtained and sequenced in whole or in part. Then, the sequence is scanned for a point at which it is desired to make a mutation (deletion, insertion or substitution) of one or more amino acids in the encoded enzyme. The sequences flanking this point are evaluated for the presence of restriction sites for replacing a short segment of the gene with an oligonucleotide pool which when expressed will encode various mutants. Such restriction sites are preferably unique sites within the protein gene so as to facilitate the replacement of the gene segment. However, any convenient restriction site that is not overly redundant in the enzyme gene may be used, provided the gene fragments generated by restriction digestion can be reassembled in proper sequence. If restriction sites are not present at locations within a convenient distance from the selected point (from 10 to 15 nucleotides), such sites are generated by substituting nucleotides in the gene in such a fashion that neither the reading frame nor the amino acids encoded are changed in the final construction. Mutation of the gene in order to change its sequence to conform to the desired sequence is accomplished by M13 primer extension in accord with generally known methods. The task of locating suitable flanking regions and evaluating the needed changes to arrive at two convenient restriction site sequences is made routine by the redundancy of the genetic code, a restriction enzyme map of the gene and the large number of different restriction enzymes. Note that if a convenient flanking restriction site is available, the above method need be used only in connection with the flanking region that does not contain a site.

[0181] Once the naturally-occurring DNA and/or synthetic DNA is cloned, the restriction sites flanking the positions to be mutated are digested with the cognate restriction enzymes and a plurality of end termini-complementary oligonucleotide cassettes are ligated into the gene. The mutagenesis is simplified by this method because all of the oligonucleotides can be synthesized so as to have the same restriction sites, and no synthetic linkers are necessary to create the restriction sites.

[0182] As used herein, "corresponding to," refers to a residue at the enumerated position in a protein or peptide, or a residue that is analogous, homologous, or equivalent to an enumerated residue in a protein or peptide. As used herein, "corresponding region," generally refers to an analogous position along related proteins or a parent protein.

[0183] As used herein, the term, "combinatorial mutagenesis" refers to methods in which libraries of variants of a starting sequence are generated. In these libraries, the variants

contain one or several mutations chosen from a predefined set of mutations. In addition, the methods provide means to introduce random mutations, which were not members of the predefined set of mutations. In some embodiments, the methods include those set forth in U.S. Application No. 09/699,250, hereby incorporated by reference. In alternative embodiments, combinatorial mutagenesis methods encompass commercially available kits (*e.g.*, QUIKCHANGE Multisite mutagenesis kit, Stratagene, San Diego, CA).

[0184] As used herein, the term “library of mutants” refers to a population of cells which are identical in most of their genome but include different homologues of one or more genes. Such libraries can be used, for example, to identify genes or operons with improved traits.

[0185] As used herein, the terms “starting gene” and “parent gene” refer to a gene of interest that encodes a protein of interest that is to be improved and/or changed using the present invention.

[0186] As used herein, the terms “multiple sequence alignment” and “MSA” refer to the sequences of multiple homologs of a starting gene that are aligned using an algorithm (*e.g.*, Clustal W).

[0187] As used herein, the terms “consensus sequence” and “canonical sequence” refer to an archetypical amino acid sequence against which all variants of a particular protein or sequence of interest are compared. The terms also refer to a sequence that sets forth the nucleotides that are most often present in a DNA sequence of interest. For each position of a gene, the consensus sequence gives the amino acid that is most abundant in that position in the MSA.

[0188] As used herein, the term “consensus mutation” refers to a difference in the sequence of a starting gene and a consensus sequence. Consensus mutations are identified by comparing the sequences of the starting gene and the consensus sequence obtained from a MSA. In some embodiments, consensus mutations are introduced into the starting gene such that it becomes more similar to the consensus sequence. Consensus mutations also include amino acid changes that change an amino acid in a starting gene to an amino acid that is more frequently found in an MSA at that position relative to the frequency of that amino acid in the starting gene. Thus, the term consensus mutation comprises all single amino acid changes that replace an amino acid of the starting gene with an amino acid that is more abundant than the amino acid in the MSA.

[0189] The terms “modified sequence” and “modified genes” are used interchangeably herein to refer to a sequence that includes a deletion, insertion or interruption of naturally occurring nucleic acid sequence. In some preferred embodiments, the expression product of

the modified sequence is a truncated protein (*e.g.*, if the modification is a deletion or interruption of the sequence). In some particularly preferred embodiments, the truncated protein retains biological activity. In alternative embodiments, the expression product of the modified sequence is an elongated protein (*e.g.*, modifications comprising an insertion into the nucleic acid sequence). In some embodiments, an insertion leads to a truncated protein (*e.g.*, when the insertion results in the formation of a stop codon). Thus, an insertion may result in either a truncated protein or an elongated protein as an expression product.

[0190] As used herein, the terms “mutant sequence” and “mutant gene” are used interchangeably and refer to a sequence that has an alteration in at least one codon occurring in a host cell’s wild-type sequence. The expression product of the mutant sequence is a protein with an altered amino acid sequence relative to the wild-type. The expression product may have an altered functional capacity (*e.g.*, enhanced enzymatic activity).

[0191] The terms “mutagenic primer” or “mutagenic oligonucleotide” (used interchangeably herein) are intended to refer to oligonucleotide compositions which correspond to a portion of the template sequence and which are capable of hybridizing thereto. With respect to mutagenic primers, the primer will not precisely match the template nucleic acid, the mismatch or mismatches in the primer being used to introduce the desired mutation into the nucleic acid library. As used herein, “non-mutagenic primer” or “non-mutagenic oligonucleotide” refers to oligonucleotide compositions that match precisely to the template nucleic acid. In one embodiment of the invention, only mutagenic primers are used. In another preferred embodiment of the invention, the primers are designed so that for at least one region at which a mutagenic primer has been included, there is also non-mutagenic primer included in the oligonucleotide mixture. By adding a mixture of mutagenic primers and non-mutagenic primers corresponding to at least one of the mutagenic primers, it is possible to produce a resulting nucleic acid library in which a variety of combinatorial mutational patterns are presented. For example, if it is desired that some of the members of the mutant nucleic acid library retain their parent sequence at certain positions while other members are mutant at such sites, the non-mutagenic primers provide the ability to obtain a specific level of non-mutant members within the nucleic acid library for a given residue. The methods of the invention employ mutagenic and non-mutagenic oligonucleotides which are generally between 10-50 bases in length, more preferably about 15-45 bases in length. However, it may be necessary to use primers that are either shorter than 10 bases or longer than 50 bases to obtain the mutagenesis result desired. With respect to corresponding mutagenic and non-mutagenic primers, it is not necessary that the corresponding

oligonucleotides be of identical length, but only that there is overlap in the region corresponding to the mutation to be added.

[0192] Primers may be added in a pre-defined ratio according to the present invention. For example, if it is desired that the resulting library have a significant level of a certain specific mutation and a lesser amount of a different mutation at the same or different site, by adjusting the amount of primer added, it is possible to produce the desired biased library. Alternatively, by adding lesser or greater amounts of non-mutagenic primers, it is possible to adjust the frequency with which the corresponding mutation(s) are produced in the mutant nucleic acid library.

[0193] The terms “wild-type sequence” or “wild-type gene” are used interchangeably herein, to refer to a sequence that is native or naturally occurring in a host cell. In some embodiments, the wild-type sequence refers to a sequence of interest that is the starting point of a protein-engineering project. The wild-type sequence may encode either a homologous or heterologous protein. A homologous protein is one the host cell would produce without intervention. A heterologous protein is one that the host cell would not produce but for the intervention.

[0194] As used herein the term “lysate” refers to a solution containing the contents of lysed cells. In some embodiments, the lysate is a bacterial cell lysate (e.g., *E. coli* cells lysed using READYLYSE™ lysozyme solution from Epicentre; or *E. coli* cells lysed using a French Pressure cell).

[0195] As used herein the term “lysozyme” refers to a glycosidase that hydrolyzes the bond between N-acetyl muramic acid and N-acetyl glucosamine, thus cleaving an important polymer in the cell wall of many bacteria. Suitable lysozymes for use with the present invention include but are not limited to hen egg white lysozyme (Sigma), T4 lysozyme, recombinant non-mammalian, non-avian lysozyme (READYLYSE™), or a fungal lysozyme.

[0196] As used herein, the term “headspace” refers to the vapor/air mixture trapped above a solid or liquid sample in a sealed vessel.

[0197] As used herein, the terms “high throughput screening” and “HTS” refer to measuring isoprene in at least 96 samples in 4 hours or less. In preferred embodiments, the sample volume is less than 2 mL.

[0198] Unless otherwise noted, all component or composition levels are in reference to the active level of that component or composition, and are exclusive of impurities, for example, residual solvents or by-products, which may be present in commercially available sources.

[0199] Enzyme components weights are based on total active protein. All percentages and ratios are calculated by weight unless otherwise indicated. All percentages and ratios are calculated based on the total composition unless otherwise.

DETAILED DESCRIPTION OF THE INVENTION

[0200] Isoprene monomer is employed in the manufacture of polyisoprene and various copolymers (with isobutylene, butadiene, styrene, or other monomers). To build a strain (prokaryotic or eukaryotic) capable of producing commercially viable levels of isoprene requires optimization of the entire pathway, either MVA to isoprene or DXP to isoprene. A key enzyme in the pathway is isoprene synthase (IspS), which converts the precursor DMAPP to isoprene. The only isoprene synthases (IspS) identified to date are those from plants such as poplar, English oak and kudzu vine. Although some bacteria, such as *Bacillus subtilis*, also produce isoprene, a prokaryotic IspS has yet to be identified and the native IspS activity in *Bacillus* is not sufficient for a commercial process. The plant IspS enzymes identified to date have been partially characterized in part by expression in *E. coli* and some of the kinetic parameters of these enzymes have been determined *in vitro* with purified protein. However, the kinetic parameters (K_m , rate etc) of the native IspS enzymes are insufficient for commercial production of isoprene in a biological host.

[0201] To solve this problem as described herein, a plant IspS is expressed in a bacterial host. In addition the IspS is engineered for a change in a property of interest. Characterization of wild-type and mutant IspS is accomplished via any means or “test” suitable and is preferably based on the assessment of properties of interest. Properties of interest include but are not limited to: pH optima, temperature stability (e.g., T_m value), intracellular and extracellular solubility, K_m value, k_{cat} value, or specific activity, as well as sensitivity to potential inhibitors including substrate or product inhibition. Oxidative and proteolytic stability are also of interest. Furthermore, activation or inhibition due to metal ion effects and ionic strength is of interest. These properties and parameters can be assessed by the conversion of DMAPP to isoprene *in vitro* with purified or partially purified isoprene synthase or *in vivo* in the context of a host organism such as *E. coli* expressing the DXP pathway, the MVA pathway, or both. It is contemplated that enzymes having various degrees of stability, solubility, activity, and/or expression level in one or more of test conditions will find use in the present invention for the production of isoprene in a diversity of hosts. High throughput methods such as those described in Example 10 are required to investigate these properties in an economical manner.

[0202] The invention features compositions and methods for the production of increased amounts of isoprene. In particular, these compositions and methods increase the rate of isoprene production and increase the total amount of isoprene that is produced. The biosynthetic processes for isoprene production described herein are a desirable alternative to using natural rubber. As discussed further below, the amount of isoprene produced by cells can be greatly increased by introducing a heterologous nucleic acid encoding an isoprene synthase (IspS) polypeptide into the cells. Isoprene synthase polypeptides convert dimethylallyl diphosphate (DMAPP) into isoprene. As shown in the examples, a heterologous *Pueraria montana* (kudzu) isoprene synthase polypeptide and variants thereof was expressed in Gram-negative bacterial cells (e.g., *Escherichia coli*). Also shown in the examples and contemplated within the scope of the invention are poplar isoprene synthase polypeptide and variants thereof was expressed in Gram-negative bacterial cells (e.g., *Escherichia coli*).

[0203] Heterologous expression of a plant IspS in bacterial host cells resulted in the production of more isoprene than the corresponding cells lacking the plant IspS.

[0204] It has been shown that mutating amino-acid residues on the surface of protease enzymes can improve their activity, expression, and stability (WO2008/153925, WO2008/153934, WO2008/153935). Surprisingly, we have found that mutating amino-acid residues on the surface of a completely different enzyme, isoprene synthase, can enhance its expression, solubility, and activity. L70R is an example of such a beneficial surface mutation.

[0205] Elucidation of the three-dimensional structure of an enzyme is essential for accurately identifying amino-acid residues on its surface. Homology modeling using structures with sequences approximately 40% identical to isoprene synthase (e.g., bornyl synthase and limonene synthase, the enzymes of known structure with closest identity to isoprene synthase) can reveal gross aspects of the modeled enzyme structure, but is insufficient to precisely identify surface-exposed residues and quantify their degree of surface exposure. Surface exposure of an amino-acid residue is quantified by the percentage of solvent-accessible surface area of its side chain.

[0206] The following classes of mutations in isoprene synthase may improve solubility of the enzyme by targeting amino-acid residues that are >50% solvent-exposed, preferably >65% solvent-exposed, and most preferably >85% solvent-exposed:

Hydrophobic → positively charged, and vice versa

Hydrophobic → negatively charged, and vice versa

Hydrophobic → neutral polar, and vice versa

Neutral polar → positively charged, and vice versa

Neutral polar → negatively charged, and vice versa

Positively charged → negatively charged, and vice versa

[0207] Additionally isoprene production by cells containing a heterologous isoprene synthase nucleic acid can be enhanced by increasing the amount of a 1-deoxy-D-xylulose-5-phosphate synthase (DXS) polypeptide, and/or an isopentenyl diphosphate isomerase (IDI) polypeptide, expressed by the cells. For example, a DXS nucleic acid and/or an IDI nucleic acid can be introduced into the cells. The DXS nucleic acid may be a heterologous nucleic acid or a duplicate copy of an endogenous nucleic acid. Similarly, the IDI nucleic acid may be a heterologous nucleic acid or a duplicate copy of an endogenous nucleic acid. In some embodiments, the amount of DXS and/or IDI polypeptide is increased by replacing the endogenous DXS and/or IDI promoters or regulatory regions with other promoters and/or regulatory regions that result in greater transcription of the DXS and/or IDI nucleic acids. In some embodiments, the cells contain both a heterologous nucleic acid encoding an isoprene synthase polypeptide (*e.g.*, a plant isoprene synthase nucleic acid) and a duplicate copy of an endogenous nucleic acid encoding an isoprene synthase polypeptide.

[0208] The encoded DXS and IDI polypeptides are part of the DXP pathway for the biosynthesis of isoprene (Figure 15). DXS polypeptides convert pyruvate and D-glyceraldehyde-3-phosphate into 1-deoxy-D-xylulose-5-phosphate. While not intending to be bound by any particular theory, it is believed that increasing the amount of DXS polypeptide increases the flow of carbon through the DXP pathway, leading to greater isoprene production. IDI polypeptides catalyze the interconversion of isopentenyl diphosphate (IPP) and dimethylallyl diphosphate (DMAPP). While not intending to be bound by any particular theory, it is believed that increasing the amount of IDI polypeptide in cells increases the amount of IPP that is converted into DMAPP, which in turn is converted into isoprene.

[0209] In some embodiments the production of isoprene by cells containing a heterologous isoprene synthase nucleic acid can be augmented by increasing expression of a MVA polypeptide in the cells (Figure 15). Exemplary MVA pathway polypeptides include any of the following polypeptides: acetyl-CoA acetyltransferase (AA-CoA thiolase) polypeptides, 3-hydroxy-3-methylglutaryl-CoA synthase (HMG-CoA synthase) polypeptides, 3-hydroxy-3-methylglutaryl-CoA reductase (HMG-CoA reductase) polypeptides, mevalonate kinase (MVK) polypeptides, phosphomevalonate kinase (PMK) polypeptides, diphosphomevalonate decarboxylase (MVD) polypeptides, IDI polypeptides, and polypeptides (*e.g.*, fusion polypeptides) having an activity of two or more MVA pathway polypeptides. For example,

one or more MVA pathway nucleic acids can be introduced into the cells. In some embodiments, the cells contain the upper MVA pathway, which includes AA-CoA thiolase, HMG-CoA synthase, and HMG-CoA reductase nucleic acids. In some embodiments, the cells contain the lower MVA pathway, which includes MVK, PMK, MVD, and IDI nucleic acids. In some embodiments, the cells contain the entire MVA pathway, which includes AA-CoA thiolase, HMG-CoA synthase, HMG-CoA reductase, MVK, PMK, MVD, and IDI nucleic acids. The MVA pathway nucleic acids may be heterologous nucleic acids or duplicate copies of endogenous nucleic acids. In some embodiments, the amount of one or more MVA pathway polypeptides is increased by replacing the endogenous promoters or regulatory regions for the MVA pathway nucleic acids with other promoters and/or regulatory regions that result in greater transcription of the MVA pathway nucleic acids. In some embodiments, the cells contain both a heterologous nucleic acid encoding an isoprene synthase polypeptide (*e.g.*, a plant isoprene synthase nucleic acid) and a duplicate copy of an endogenous nucleic acid encoding an isoprene synthase polypeptide.

[0210] In some embodiments, at least a portion of the cells maintain the heterologous isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acid for at least about 5, 10, 20, 50, 75, 100, 200, 300, or more cell divisions in a continuous culture (such as a continuous culture without dilution). In some embodiments of any of the aspects of the invention, the nucleic acid comprising the heterologous or duplicate copy of an endogenous isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acid also comprises a selective marker, such as a kanamycin, ampicillin, carbenicillin, gentamicin, hygromycin, phleomycin, bleomycin, neomycin, or chloramphenicol antibiotic resistance nucleic acid.

I. Exemplary Polypeptides and Nucleic Acids

[0211] Various isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides and nucleic acids can be used in the compositions and methods of the invention.

[0212] As used herein, “polypeptides” includes polypeptides, proteins, peptides, fragments of polypeptides, and fusion polypeptides that include part or all of a first polypeptide (*e.g.*, an isoprene synthase, DXS, IDI, or MVA pathway polypeptide) and part or all of a second polypeptide (*e.g.*, a peptide that facilitates purification or detection of the fusion polypeptide, such as a His-tag). In some embodiments, the fusion polypeptide has an activity of two or more MVA pathway polypeptides (such as AA-CoA thiolase and HMG-CoA reductase polypeptides). In some embodiments, the polypeptide is a naturally-occurring polypeptide

(such as the polypeptide encoded by an *Enterococcus faecalis mvaE* nucleic acid) that has an activity of two or more MVA pathway polypeptides.

[0213] In various embodiments, a polypeptide has at least or about 50, 100, 150, 175, 200, 250, 300, 350, 400, or more amino acids. In some embodiments, the polypeptide fragment contains at least or about 25, 50, 75, 100, 150, 200, 300, or more contiguous amino acids from a full-length polypeptide and has at least or about 5%, 10%, 20%, 30%, 40%, 50%, 60%, 70%, 80%, 90%, 95%, or 100% of an activity of a corresponding full-length polypeptide. In particular embodiments the polypeptide includes a segment of or the entire amino acid sequence of any naturally-occurring isoprene synthase, DXS, IDI, or MVA pathway polypeptide. In some embodiments, the polypeptide has one or more mutations compared to the sequence of a wild-type (*i.e.*, a sequence occurring in nature) isoprene synthase, DXS, IDI, or MVA pathway polypeptide.

[0214] In some embodiments, the polypeptide is an isolated polypeptide. As used herein, an “isolated polypeptide” is not part of a library of polypeptides, such as a library of 2, 5, 10, 20, 50 or more different polypeptides and is separated from at least one component with which it occurs in nature. An isolated polypeptide can be obtained, for example, by expression of a recombinant nucleic acid encoding the polypeptide.

[0215] In some embodiments, the polypeptide is a heterologous polypeptide. By “heterologous polypeptide” is meant a polypeptide whose amino acid sequence is not identical to that of another polypeptide naturally expressed in the same host cell.

[0216] As used herein, a “nucleic acid” refers to two or more deoxyribonucleotides and/or ribonucleotides in either single or double-stranded form. In some embodiments, the nucleic acid is a recombinant nucleic acid. By “recombinant nucleic acid” means a nucleic acid of interest that is free of one or more nucleic acids (*e.g.*, genes), which in the genome occurring in nature of the organism from which the nucleic acid of interest is derived, flank the nucleic acid of interest. The term therefore includes, for example, a recombinant DNA which is incorporated into a vector, into an autonomously replicating plasmid or virus, or into the genomic DNA of a prokaryote or eukaryote, or which exists as a separate molecule (*e.g.*, a cDNA, a genomic DNA fragment, or a cDNA fragment produced by PCR or restriction endonuclease digestion) independent of other sequences.

[0217] In various embodiments, the nucleic acid is a recombinant nucleic acid. For instance, in some embodiments, an isoprene synthase, DXS, IDI, or MVA pathway nucleic acid is operably linked to another nucleic acid encoding all or a portion of another polypeptide such that the recombinant nucleic acid encodes a fusion polypeptide that includes

an isoprene synthase, DXS, IDI, or MVA pathway polypeptide and all or part of another polypeptide (*e.g.*, a peptide that facilitates purification or detection of the fusion polypeptide, such as a His-tag). In some embodiments, part or all of a recombinant nucleic acid is chemically synthesized. In some embodiments, the nucleic acid is a heterologous nucleic acid. By “heterologous nucleic acid” is meant a nucleic acid whose nucleic acid sequence is not identical to that of another nucleic acid naturally found in the same host cell.

[0218] In particular embodiments the nucleic acid includes a segment of or the entire nucleic acid sequence of any naturally-occurring isoprene synthase, DXS, IDI, or MVA pathway nucleic acid. In some embodiments, the nucleic acid includes at least or about 50, 100, 150, 200, 300, 400, 500, 600, 700, 800, or more contiguous nucleotides from a naturally-occurring isoprene synthase nucleic acid DXS, IDI, or MVA pathway nucleic acid. In some embodiments, the nucleic acid has one or more mutations compared to the sequence of a wild-type (*i.e.*, a sequence occurring in nature) isoprene synthase, DXS, IDI, or MVA pathway nucleic acid. In some embodiments, the nucleic acid has one or more mutations (*e.g.*, a silent mutation) that increase the transcription or translation of isoprene synthase, DXS, IDI, or MVA pathway nucleic acid. In some embodiments, the nucleic acid is a degenerate variant of any nucleic acid encoding an isoprene synthase, DXS, IDI, or MVA pathway polypeptide.

[0219] “Codon degeneracy” refers to divergence in the genetic code permitting variation of the nucleotide sequence without affecting the amino acid sequence of an encoded polypeptide. The skilled artisan is well aware of the “codon-bias” exhibited by a specific host cell in usage of nucleotide codons to specify a given amino acid. Therefore, when synthesizing a nucleic acid for improved expression in a host cell, it is desirable in some embodiments to design the nucleic acid such that its frequency of codon usage approaches the frequency of preferred codon usage of the host cell.

[0220] The accession numbers of exemplary isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides and nucleic acids are listed in Appendix 1 of U.S. Application No. 61/013,574, herein incorporated by reference in its entirety, particularly with respect to the amino acid and nucleic acid sequences of isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides and nucleic acids). The Kegg database also contains the amino acid and nucleic acid sequences of numerous exemplary isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides and nucleic acids (*See, e.g.*, the world-wide web at “genome.jp/kegg/pathway/map/map00100.html” and the sequences therein, which are each hereby incorporated by reference in their entireties, particularly with respect to the amino

acid and nucleic acid sequences of isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides and nucleic acids). In some embodiments, one or more of the isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides and/or nucleic acids have a sequence identical to a sequence publicly available on December 12, 2007, such as any of the sequences that correspond to any of the accession numbers in Appendix 1 of U.S. Application No. 61/013,574, or any of the sequences present in the Kegg database as of the date of this filing. Additional exemplary isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides and nucleic acids are described further below.

Exemplary Isoprene Synthase Polypeptides and Nucleic Acids

[0221] As noted above, isoprene synthase polypeptides convert dimethylallyl diphosphate (DMAPP) into isoprene. Exemplary isoprene synthase polypeptides include polypeptides, fragments of polypeptides, peptides, and fusions polypeptides that have at least one activity of an isoprene synthase polypeptide. Standard methods can be used to determine whether a polypeptide has isoprene synthase polypeptide activity by measuring the ability of the polypeptide to convert DMAPP into isoprene *in vitro*, in a cell extract, or *in vivo*. Isoprene synthase polypeptide activity in the cell extract can be measured, for example, as described in Silver *et al.*, J. Biol. Chem. 270:13010-13016, 1995 and references therein, which are each hereby incorporated by reference in their entireties, particularly with respect to assays for isoprene synthase polypeptide activity. DMAPP (Sigma) is evaporated to dryness under a stream of nitrogen and rehydrated to a concentration of 100 mM in 100 mM potassium phosphate buffer pH 8.2 and stored at -20 °C. To perform the assay, a solution of 5 µl of 1M MgCl₂, 1 mM (250 µg/ml) DMAPP, 65 µl of Plant Extract Buffer (PEB) (50 mM Tris-HCl, pH 8.0, 20 mM MgCl₂, 5% glycerol, and 2 mM DTT) is added to 25 µl of cell extract in a 20 ml Headspace vial with a metal screw cap and teflon coated silicon septum (Agilent Technologies) and cultured at 37 °C for 15 minutes with shaking. The reaction is quenched by adding 200 µl of 250 mM EDTA or by heat inactivation, and isoprene is quantified by GC/MS.

[0222] Exemplary isoprene synthase nucleic acids include nucleic acids that encode a polypeptide, fragment of a polypeptide, peptide, or fusion polypeptide that has at least one activity of an isoprene synthase polypeptide. Exemplary isoprene synthase polypeptides and nucleic acids include naturally-occurring polypeptides and nucleic acids from any of the

source organisms described herein as well as mutant polypeptides and nucleic acids derived from any of the source organisms described herein.

[0223] In some embodiments, the isoprene synthase polypeptide or nucleic acid is from the family Fabaceae, the family Salicaceae, or the family Fagaceae. In some embodiments, the isoprene synthase polypeptide or nucleic acid is a naturally-occurring polypeptide or nucleic acid from *Pueraria montana* (kudzu) (Sharkey *et al.*, Plant Physiology 137: 700-712, 2005), poplar (such as *Populus alba x tremula* CAC35696) Miller *et al.*, Planta 213: 483-487, 2001) aspen (such as *Populus tremuloides*) Silver *et al.*, JBC 270(22): 13010-1316, 1995), or English Oak (*Quercus robur*) (Zimmer *et al.*, WO 98/02550), which are each hereby incorporated by reference in their entireties, particularly with respect to isoprene synthase nucleic acids and the expression of isoprene synthase polypeptides. Suitable isoprene synthases include, but are not limited to, those identified by GenBank Accession Nos. AY341431, AY316691, AY279379, AJ457070, and AY182241, which are each hereby incorporated by reference in their entireties, particularly with respect to sequences of isoprene synthase nucleic acids and polypeptides. In some embodiments, the isoprene synthase polypeptide or nucleic acid is not a naturally-occurring polypeptide or nucleic acid from *Quercus robur* (*i.e.*, the isoprene synthase polypeptide or nucleic acid is an isoprene synthase polypeptide or nucleic acid other than a naturally-occurring polypeptide or nucleic acid from *Quercus robur*). In some embodiments, the isoprene synthase nucleic acid or polypeptide is not a naturally-occurring polypeptide or nucleic acid from poplar (such as *Populus alba x tremula* CAC35696).

Exemplary DXS Polypeptides and Nucleic Acids

[0224] As noted above, 1-deoxy-D-xylulose-5-phosphate synthase (DXS) polypeptides convert pyruvate and D-glyceraldehyde-3-phosphate into 1-deoxy-D-xylulose-5-phosphate. Exemplary DXS polypeptides include polypeptides, fragments of polypeptides, peptides, and fusions polypeptides that have at least one activity of a DXS polypeptide. Standard methods (such as those described herein) can be used to determine whether a polypeptide has DXS polypeptide activity by measuring the ability of the polypeptide to convert pyruvate and D-glyceraldehyde-3-phosphate into 1-deoxy-D-xylulose-5-phosphate *in vitro*, in a cell extract, or *in vivo*. Exemplary DXS nucleic acids include nucleic acids that encode a polypeptide, fragment of a polypeptide, peptide, or fusion polypeptide that has at least one activity of a DXS polypeptide. Exemplary DXS polypeptides and nucleic acids include naturally-occurring polypeptides and nucleic acids from any of the source organisms described herein

as well as mutant polypeptides and nucleic acids derived from any of the source organisms described herein.

Exemplary IDI Polypeptides and Nucleic Acids

[0225] Isopentenyl diphosphate isomerase polypeptides (isopentenyl-diphosphate delta-isomerase or IDI) catalyzes the interconversion of isopentenyl diphosphate (IPP) and dimethylallyl diphosphate (DMAPP) (*e.g.*, converting IPP into DMAPP and/or converting DMAPP into IPP). Exemplary IDI polypeptides include polypeptides, fragments of polypeptides, peptides, and fusions polypeptides that have at least one activity of an IDI polypeptide. Standard methods (such as those described herein) can be used to determine whether a polypeptide has IDI polypeptide activity by measuring the ability of the polypeptide to interconvert IPP and DMAPP *in vitro*, in a cell extract, or *in vivo*. Exemplary IDI nucleic acids include nucleic acids that encode a polypeptide, fragment of a polypeptide, peptide, or fusion polypeptide that has at least one activity of an IDI polypeptide. Exemplary IDI polypeptides and nucleic acids include naturally-occurring polypeptides and nucleic acids from any of the source organisms described herein as well as mutant polypeptides and nucleic acids derived from any of the source organisms described herein.

Exemplary MVA Pathway Polypeptides and Nucleic Acids

[0226] Exemplary MVA pathway polypeptides include acetyl-CoA acetyltransferase (AA-CoA thiolase) polypeptides, 3-hydroxy-3-methylglutaryl-CoA synthase (HMG-CoA synthase) polypeptides, 3-hydroxy-3-methylglutaryl-CoA reductase (HMG-CoA reductase) polypeptides, mevalonate kinase (MVK) polypeptides, phosphomevalonate kinase (PMK) polypeptides, diphosphomevalonate decarboxylase (MVD) polypeptides, IDI polypeptides, and polypeptides (*e.g.*, fusion polypeptides) having an activity of two or more MVA pathway polypeptides. In particular, MVA pathway polypeptides include polypeptides, fragments of polypeptides, peptides, and fusions polypeptides that have at least one activity of an MVA pathway polypeptide. Exemplary MVA pathway nucleic acids include nucleic acids that encode a polypeptide, fragment of a polypeptide, peptide, or fusion polypeptide that has at least one activity of an MVA pathway polypeptide. Exemplary MVA pathway polypeptides and nucleic acids include naturally-occurring polypeptides and nucleic acids from any of the source organisms described herein as well as mutant polypeptides and nucleic acids derived from any of the source organisms described herein.

[0227] In particular, acetyl-CoA acetyltransferase polypeptides (AA-CoA thiolase or AACT) convert two molecules of acetyl-CoA into acetoacetyl-CoA. Standard methods (such as those described herein) can be used to determine whether a polypeptide has AA-CoA thiolase polypeptide activity by measuring the ability of the polypeptide to convert two molecules of acetyl-CoA into acetoacetyl-CoA *in vitro*, in a cell extract, or *in vivo*.

[0228] 3-hydroxy-3-methylglutaryl-CoA synthase (HMG-CoA synthase or HMGS) polypeptides convert acetoacetyl-CoA into 3-hydroxy-3-methylglutaryl-CoA. Standard methods (such as those described herein) can be used to determine whether a polypeptide has HMG-CoA synthase polypeptide activity by measuring the ability of the polypeptide to convert acetoacetyl-CoA into 3-hydroxy-3-methylglutaryl-CoA *in vitro*, in a cell extract, or *in vivo*.

[0229] 3-hydroxy-3-methylglutaryl-CoA reductase (HMG-CoA reductase or HMGR) polypeptides convert 3-hydroxy-3-methylglutaryl-CoA into mevalonate. Standard methods (such as those described herein) can be used to determine whether a polypeptide has HMG-CoA reductase polypeptide activity by measuring the ability of the polypeptide to convert 3-hydroxy-3-methylglutaryl-CoA into mevalonate *in vitro*, in a cell extract, or *in vivo*.

[0230] Mevalonate kinase (MVK) polypeptides phosphorylates mevalonate to form mevalonate-5-phosphate. Standard methods (such as those described herein) can be used to determine whether a polypeptide has MVK polypeptide activity by measuring the ability of the polypeptide to convert mevalonate into mevalonate-5-phosphate *in vitro*, in a cell extract, or *in vivo*.

[0231] Phosphomevalonate kinase (PMK) polypeptides phosphorylates mevalonate-5-phosphate to form mevalonate-5-diphosphate. Standard methods (such as those described herein) can be used to determine whether a polypeptide has PMK polypeptide activity by measuring the ability of the polypeptide to convert mevalonate-5-phosphate into mevalonate-5-diphosphate *in vitro*, in a cell extract, or *in vivo*.

[0232] Diphosphomevalonate decarboxylase (MVD or DPMDC) polypeptides convert mevalonate-5-diphosphate into isopentenyl diphosphate polypeptides (IPP). Standard methods (such as those described herein) can be used to determine whether a polypeptide has MVD polypeptide activity by measuring the ability of the polypeptide to convert mevalonate-5-diphosphate into IPP *in vitro*, in a cell extract, or *in vivo*.

Exemplary Methods for Isolating Nucleic Acids

[0233] Isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acids can be isolated using standard methods. Methods of obtaining desired nucleic acids from a source organism of interest (such as a bacterial genome) are common and well known in the art of molecular biology (*see*, for example, WO 2004/033646 and references cited therein, which are each hereby incorporated by reference in their entireties, particularly with respect to the isolation of nucleic acids of interest). For example, if the sequence of the nucleic acid is known (such as any of the known nucleic acids described herein), suitable genomic libraries may be created by restriction endonuclease digestion and may be screened with probes complementary to the desired nucleic acid sequence. Once the sequence is isolated, the DNA may be amplified using standard primer directed amplification methods such as polymerase chain reaction (PCR) (U.S. Patent No. 4,683,202, which is incorporated by reference in its entirety, particularly with respect to PCR methods) to obtain amounts of DNA suitable for transformation using appropriate vectors. Alternatively, isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acids (such as any isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acids with a known nucleic acid sequence) can be chemically synthesized using standard methods.

[0234] Additional isoprene synthase, DXS, IDI, or MVA pathway polypeptides and nucleic acids that may be suitable for use in the compositions and methods described herein can be identified using standard methods. For example, cosmid libraries of the chromosomal DNA of organisms known to produce isoprene naturally can be constructed in organisms such as *E. coli*, and then screened for isoprene production. Additional methods for obtaining isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acids include screening a metagenomic library by assay (such as the headspace assay described herein) or by PCR using primers directed against nucleotides encoding for a length of conserved amino acids (for example, at least 3 conserved amino acids). Conserved amino acids can be identified by aligning amino acid sequences of known isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides. Conserved amino acids for isoprene synthase polypeptides can be identified based on aligned sequences of known isoprene synthase polypeptides. An organism found to produce isoprene naturally can be subjected to standard protein purification methods (which are well known in the art) and the resulting purified polypeptide can be sequenced using standard methods. Other methods are found in the literature (See, e.g., Julsing et al., *Applied. Microbiol. Biotechnol.* 75: 1377-84, 2007; and Withers et al., *Appl Environ Microbiol.* 73:6277-83, 2007, which are each hereby incorporated by reference in their entireties, particularly with respect to identification of nucleic acids involved in the synthesis of isoprene).

[0235] Additionally, standard sequence alignment and/or structure prediction programs can be used to identify additional DXS, IDI, or MVA pathway polypeptides and nucleic acids based on the similarity of their primary and/or predicted polypeptide secondary structure with that of known DXS, IDI, or MVA pathway polypeptides and nucleic acids. Standard databases such as the swissprot-trembl database (world-wide web at “expasy.org”, Swiss Institute of Bioinformatics Swiss-Prot group CMU - 1 rue Michel Servet CH-1211 Geneva 4, Switzerland) can also be used to identify isoprene synthase, DXS, IDI, or MVA pathway polypeptides and nucleic acids. The secondary and/or tertiary structure of an isoprene synthase, DXS, IDI, or MVA pathway polypeptide can be predicted using the default settings of standard structure prediction programs, such as PredictProtein (Rost et al., The PredictProtein Server. Nucleic Acids Research 32(Web Server issue):W321-W326, 2004). Alternatively, the actual secondary and/or tertiary structure of an isoprene synthase, DXS, IDI, or MVA pathway polypeptide can be determined using standard methods. Additional isoprene synthase, DXS, IDI, or MVA pathway nucleic acids can also be identified by hybridization to probes generated from known isoprene synthase, DXS, IDI, or MVA pathway nucleic acids.

Exemplary Promoters and Vectors

[0236] Any of the isoprene synthase, DXS, IDI, or MVA pathway nucleic acid described herein can be included in one or more vectors. Accordingly, the invention also features vectors with one more nucleic acids encoding any of the isoprene synthase, DXS, IDI, or MVA pathway polypeptides that are described herein. As used herein, a “vector” means a construct that is capable of delivering, and desirably expressing one or more nucleic acids of interest in a host cell. Examples of vectors include, but are not limited to, plasmids, viral vectors, DNA or RNA expression vectors, cosmids, and phage vectors. In some embodiments, the vector contains a nucleic acid under the control of an expression control sequence.

[0237] As used herein, an “expression control sequence” means a nucleic acid sequence that directs transcription of a nucleic acid of interest. An expression control sequence can be a promoter, such as a constitutive or an inducible promoter, or an enhancer. An “inducible promoter” is a promoter that is active under environmental or developmental regulation. The expression control sequence is operably linked to the nucleic acid segment to be transcribed.

[0238] In some embodiments, the vector contains a selective marker. The term “selective marker” refers to a nucleic acid capable of expression in a host cell that allows for ease of

selection of those host cells containing an introduced nucleic acid or vector. Examples of selectable markers include, but are not limited to, antibiotic resistance nucleic acids (*e.g.*, kanamycin, ampicillin, carbenicillin, gentamicin, hygromycin, phleomycin, bleomycin, neomycin, or chloramphenicol) and/or nucleic acids that confer a metabolic advantage, such as a nutritional advantage on the host cell. In some embodiments, an isoprene synthase, DXS, IDI, or MVA pathway nucleic acid integrates into a chromosome of the cells without a selective marker.

[0239] Suitable vectors are those that are compatible with the host cell employed. Suitable vectors can be derived, for example, from a bacterium, a virus (such as bacteriophage T7 or a M-13 derived phage), a cosmid, a yeast, or a plant. Protocols for obtaining and using such vectors are known in the art (*See, e.g.*, Sambrook *et al.*, Molecular Cloning: A Laboratory Manual, 2nd ed., Cold Spring Harbor, 1989, which is hereby incorporated by reference in its entirety, particularly with respect to the use of vectors).

[0240] Promoters are well known in the art. Any promoter that functions in the host cell can be used for expression of an isoprene synthase, DXS, IDI, or MVA pathway nucleic acid in the host cell. Initiation control regions or promoters, which are useful to drive expression of isoprene synthase, DXS, IDI, or MVA pathway nucleic acids in various host cells are numerous and familiar to those skilled in the art (*see, for example*, WO 2004/033646 and references cited therein, which are each hereby incorporated by reference in their entireties, particularly with respect to vectors for the expression of nucleic acids of interest). Virtually any promoter capable of driving these nucleic acids is suitable for the present invention including, but not limited to *lac*, *trp*, λP_L , λP_R , T7, *tac*, and *trc* (useful for expression in *E. coli*).

[0241] In some embodiments, a glucose isomerase promoter is used (*see, for example*, U.S. Patent No. 7,132,527 and references cited therein, which are each hereby incorporated by reference in their entireties, particularly with respect promoters and plasmid systems for expressing polypeptides of interest). Reported glucose isomerase promoter mutants can be used to vary the level of expression of the polypeptide encoded by a nucleic acid operably linked to the glucose isomerase promoter (U.S. Patent No. 7,132,527). In various embodiments, the glucose isomerase promoter is contained in a low, medium, or high copy plasmid (U.S. Patent No. 7,132,527).

[0242] In various embodiments, an isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acid is contained in a low copy plasmid (*e.g.*, a plasmid that is maintained at about 1

to about 4 copies per cell), medium copy plasmid (*e.g.*, a plasmid that is maintained at about 10 to about 15 copies per cell), or high copy plasmid (*e.g.*, a plasmid that is maintained at about 50 or more copies per cell). In some embodiments, the heterologous or extra endogenous isoprene synthase, DXS, IDI, or MVA pathway nucleic acid is operably linked to a T7 promoter. In some embodiments, the heterologous or extra endogenous isoprene synthase, DXS, IDI, or MVA pathway nucleic acid operably linked to a T7 promoter is contained in a medium or high copy plasmid. In some embodiments, the heterologous or extra endogenous isoprene synthase, DXS, IDI, or MVA pathway nucleic acid is operably linked to a Trc promoter. In some embodiments, the heterologous or extra endogenous isoprene synthase, DXS, IDI, or MVA pathway nucleic acid operably linked to a Trc promoter is contained in a medium or high copy plasmid. In some embodiments, the heterologous or extra endogenous isoprene synthase, DXS, IDI, or MVA pathway nucleic acid is operably linked to a Lac promoter. In some embodiments, the heterologous or extra endogenous isoprene synthase, DXS, IDI, or MVA pathway nucleic acid operably linked to a Lac promoter is contained in a low copy plasmid. In some embodiments, the heterologous or extra endogenous isoprene synthase, DXS, IDI, or MVA pathway nucleic acid is operably linked to an endogenous promoter, such as an endogenous *Escherichia*, *Pantaea*, *Bacillus*, *Yarrowia*, *Streptomyces*, or *Trichoderma* promoter or an endogenous alkaline serine protease, isoprene synthase, DXS, IDI, or MVA pathway promoter. In some embodiments, the heterologous or extra endogenous isoprene synthase, DXS, IDI, or MVA pathway nucleic acid operably linked to an endogenous promoter is contained in a high copy plasmid. In some embodiments, the vector is a replicating plasmid that does not integrate into a chromosome in the cells. In some embodiments part or all of the vector integrates into a chromosome in the cells.

[0243] In some embodiments, the expression vector also includes a termination sequence. Termination control regions may also be derived from various genes native to the host cell. In some embodiments, the termination sequence and the promoter sequence are derived from the same source. In another embodiment, the termination sequence is endogenous to the host cell.

[0244] In some embodiments, the promoter, coding, region, and terminator all originate from the isoprene synthase, DXS, IDI, or MVA pathway nucleic acid to be expressed. In some embodiments, the coding region for an isoprene synthase, DXS, IDI, or MVA pathway nucleic acid is inserted into a general-purpose expression vector such that it is under the transcriptional control of the expression construct promoter and terminator sequences. In

some embodiments, genes or part thereof are inserted downstream of the strong *cbh1* promoter.

[0245] An isoprene synthase, DXS, IDI, or MVA pathway nucleic acid can be incorporated into a vector, such as an expression vector, using standard techniques (Sambrook *et al.*, Molecular Cloning: A Laboratory Manual, Cold Spring Harbor, 1982, which is hereby incorporated by reference in its entirety, particularly with respect to the screening of appropriate DNA sequences and the construction of vectors). Methods used to ligate the DNA construct comprising a nucleic acid of interest (such as an isoprene synthase, DXS, IDI, or MVA pathway nucleic acid), a promoter, a terminator, and other sequences and to insert them into a suitable vector are well known in the art. For example, restriction enzymes can be used to cleave the isoprene synthase, DXS, IDI, or MVA pathway nucleic acid and the vector. Then, the compatible ends of the cleaved isoprene synthase, DXS, IDI, or MVA pathway nucleic acid and the cleaved vector can be ligated. Linking is generally accomplished by ligation at convenient restriction sites. If such sites do not exist, the synthetic oligonucleotide linkers are used in accordance with conventional practice (*see*, Sambrook *et al.*, Molecular Cloning: A Laboratory Manual, 2nd ed., Cold Spring Harbor, 1989, and Bennett and Lasure, More Gene Manipulations in Fungi, Academic Press, San Diego, pp 70–76, 1991, which are each hereby incorporated by reference in their entireties, particularly with respect to oligonucleotide linkers). Additionally, vectors can be constructed using known recombination techniques (*e.g.*, Invitrogen Life Technologies, Gateway Technology).

[0246] In some embodiments, it may be desirable to over-express isoprene synthase, DXS, IDI, or MVA pathway nucleic acids at levels far higher than currently found in naturally-occurring cells. This result may be accomplished by the selective cloning of the nucleic acids encoding those polypeptides into multicopy plasmids or placing those nucleic acids under a strong inducible or constitutive promoter. Methods for over-expressing desired polypeptides are common and well known in the art of molecular biology and examples may be found in Sambrook *et al.*, Molecular Cloning: A Laboratory Manual, 2nd ed., Cold Spring Harbor, 1989, which is hereby incorporated by reference in its entirety, particularly with respect to cloning techniques.

[0247] The following resources include descriptions of additional general methodology useful in accordance with the invention: Kreigler, Gene Transfer and Expression; A Laboratory Manual, 1990; and Ausubel *et al.*, Eds. Current Protocols in Molecular Biology,

1994, which are each hereby incorporated by reference in their entireties, particularly with respect to molecular biology and cloning techniques.

Exemplary Source Organisms

[0248] Isoprene synthase, DXS, IDI, or MVA pathway nucleic acids (and their encoded polypeptides) can be obtained from any organism that naturally contains isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acids. As noted above, isoprene is formed naturally by a variety of organisms, such as bacteria, yeast, plants, and animals. Organisms contain the MVA pathway, DXP pathway, or both the MVA and DXP pathways for producing isoprene (Figure 15). Thus, DXS nucleic acids can be obtained, *e.g.*, from any organism that contains the DXP pathway or contains both the MVA and DXP pathways. IDI and isoprene synthase nucleic acids can be obtained, *e.g.*, from any organism that contains the MVA pathway, DXP pathway, or both the MVA and DXP pathways. MVA pathway nucleic acids can be obtained, *e.g.*, from any organism that contains the MVA pathway or contains both the MVA and DXP pathways.

[0249] In some embodiments, the nucleic acid sequence of the isoprene synthase, DXS, IDI, or MVA pathway nucleic acid is identical to the sequence of a nucleic acid that is produced by any of the following organisms in nature. In some embodiments, the amino acid sequence of the isoprene synthase, DXS, IDI, or MVA pathway polypeptide is identical to the sequence of a polypeptide that is produced by any of the following organisms in nature. In some embodiments, the isoprene synthase, DXS, IDI, or MVA pathway nucleic acid or polypeptide is a mutant nucleic acid or polypeptide derived from any of the organisms described herein. As used herein, “derived from” refers to the source of the nucleic acid or polypeptide into which one or more mutations is introduced. For example, a polypeptide that is “derived from a plant polypeptide” refers to polypeptide of interest that results from introducing one or more mutations into the sequence of a wild-type (*i.e.*, a sequence occurring in nature) plant polypeptide.

[0250] In some embodiments, the source organism is a bacterium, such as strains of *Escherichia* (*e.g.*, *E. coli*), or strains of *Bacillus* (*e.g.*, *B. subtilis*).

[0251] As used herein, “the genus *Escherichia*” includes all species within the genus “*Escherichia*,” as known to those of skill in the art, including but not limited to *E. coli*, *E. adecarboxylata*, *E. albertii*, *E. blattae*, *E. fergusonii*, *E. hermannii*, *E. senegalensis*, and *E. vulneris*. The genus “*Escherichia*” is defined as Gram-negative, non-spore forming,

facultatively anaerobic, rod-shaped bacteria are classified as members of the Family *Enterobacteriaceae*, Order *Enterobacteriales*, Class *Gamma Proteobacteria*.

[0252] As used herein, “the genus *Bacillus*” includes all species within the genus “*Bacillus*,” as known to those of skill in the art, including but not limited to *B. subtilis*, *B. licheniformis*, *B. lentus*, *B. brevis*, *B. stearothermophilus*, *B. alkalophilus*, *B. amyloliquefaciens*, *B. clausii*, *B. halodurans*, *B. megaterium*, *B. coagulans*, *B. circulans*, *B. lautus*, and *B. thuringiensis*. It is recognized that the genus *Bacillus* continues to undergo taxonomical reorganization. Thus, it is intended that the genus include species that have been reclassified, including but not limited to such organisms as *B. stearothermophilus*, which is now named “*Geobacillus stearothermophilus*.” The production of resistant endospores in the presence of oxygen is considered the defining feature of the genus *Bacillus*, although this characteristic also applies to the recently named *Alicyclobacillus*, *Amphibacillus*, *Aneurinibacillus*, *Anoxybacillus*, *Brevibacillus*, *Filobacillus*, *Gracilibacillus*, *Halobacillus*, *Paenibacillus*, *Salibacillus*, *Thermobacillus*, *Ureibacillus*, and *Virgibacillus*.

Exemplary Host Cells

[0253] A variety of host cells can be used to express isoprene synthase, DXS, IDI, and/or MVA pathway polypeptides and to produce isoprene in the methods of the claimed invention. Exemplary host cells include cells from any of the organisms listed in the prior section under the heading “*Exemplary Source Organisms*.” The host cell may be a cell that naturally produces isoprene or a cell that does not naturally produce isoprene. In some embodiments, the host cell naturally produces isoprene using the DXP pathway, and an isoprene synthase, DXS, and/or IDI nucleic acid is added to enhance production of isoprene using this pathway. In some embodiments, the host cell naturally produces isoprene using the MVA pathway, and an isoprene synthase and/or one or more MVA pathway nucleic acids are added to enhance production of isoprene using this pathway. In some embodiments, the host cell naturally produces isoprene using the DXP pathway and one or more MVA pathway nucleic acids are added to produce isoprene using part or all of the MVA pathway as well as the DXP pathway. In some embodiments, the host cell naturally produces isoprene using both the DXP and MVA pathways and one or more isoprene synthase, DXS, IDI, or MVA pathway nucleic acids are added to enhance production of isoprene by one or both of these pathways.

Exemplary Transformation Methods

[0254] Isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acids or vectors containing them can be inserted into a host cell (*e.g.*, a bacterial cell) using standard techniques for expression of the encoded isoprene synthase, DXS, IDI, and/or MVA pathway polypeptide. Introduction of a DNA construct or vector into a host cell can be performed using techniques such as transformation, electroporation, nuclear microinjection, transduction, transfection (*e.g.*, lipofection mediated or DEAE-Dextrin mediated transfection or transfection using a recombinant phage virus), incubation with calcium phosphate DNA precipitate, high velocity bombardment with DNA-coated microprojectiles, and protoplast fusion. General transformation techniques are known in the art (*see, e.g.*, Current Protocols in Molecular Biology (F. M. Ausubel *et al.* (eds) Chapter 9, 1987; Sambrook *et al.*, Molecular Cloning: A Laboratory Manual, 2nd ed., Cold Spring Harbor, 1989; and Campbell *et al.*, *Curr Genet*, 16:53-56, 1989, which are each hereby incorporated by reference in their entireties, particularly with respect to transformation methods). The introduced nucleic acids may be integrated into chromosomal DNA or maintained as extrachromosomal replicating sequences.

Exemplary Cell Culture Media

[0255] The invention also includes a cell or a population of cells in culture that produce isoprene. By “cells in culture” is meant two or more cells in a solution (*e.g.*, a cell medium) that allows the cells to undergo one or more cell divisions. “Cells in culture” do not include plant cells that are part of a living, multicellular plant containing cells that have differentiated into plant tissues. In various embodiments, the cell culture includes at least or about 10, 20, 50, 100, 200, 500, 1,000, 5,000, 10,000 or more cells.

[0256] Any carbon source can be used to cultivate the host cells. The term “carbon source” refers to one or more carbon-containing compounds capable of being metabolized by a host cell or organism. For example, the cell medium used to cultivate the host cells may include any carbon source suitable for maintaining the viability or growing the host cells.

[0257] In some embodiments, the carbon source is a carbohydrate (such as monosaccharide, disaccharide, oligosaccharide, or polysaccharides), invert sugar (*e.g.*, enzymatically treated sucrose syrup), glycerol, glycerine (*e.g.*, a glycerine byproduct of a biodiesel or soap-making process), dihydroxyacetone, one-carbon source, fatty acid (*e.g.*, a saturated fatty acid, unsaturated fatty acid, or polyunsaturated fatty acid), lipid, phospholipid, glycerolipid, monoglyceride, diglyceride, triglyceride, polypeptide (*e.g.*, a microbial or plant protein or peptide), renewable carbon source (*e.g.*, a biomass carbon source such as a hydrolyzed biomass carbon source; beet sugar or cane sugar molasses), yeast extract,

component from a yeast extract, polymer, acid, alcohol, aldehyde, ketone, amino acid, succinate, lactate, acetate, ethanol, or any combination of two or more of the foregoing. In some embodiments, the carbon source is a product of photosynthesis, including, but not limited to, glucose.

[0258] Exemplary monosaccharides include glucose and fructose; exemplary oligosaccharides include lactose and sucrose, and exemplary polysaccharides include starch and cellulose. Exemplary carbohydrates include C6 sugars (*e.g.*, fructose, mannose, galactose, or glucose) and C5 sugars (*e.g.*, xylose or arabinose). In some embodiments, the cell medium includes a carbohydrate as well as a carbon source other than a carbohydrate (*e.g.*, glycerol, glycerine, dihydroxyacetone, one-carbon source, fatty acid, lipid, phospholipid, glycerolipid, monoglyceride, diglyceride, triglyceride, renewable carbon source, or a component from a yeast extract). In some embodiments, the cell medium includes a carbohydrate as well as a polypeptide (*e.g.*, a microbial or plant protein or peptide). In some embodiments, the microbial polypeptide is a polypeptide from yeast or bacteria. In some embodiments, the plant polypeptide is a polypeptide from soy, corn, canola, jatropha, palm, peanut, sunflower, coconut, mustard, rapeseed, cottonseed, palm kernel, olive, safflower, sesame, or linseed.

[0259] In some embodiments, the concentration of the carbohydrate is at least or about 5 grams per liter of broth (g/L, wherein the volume of broth includes both the volume of the cell medium and the volume of the cells), such as at least or about 10, 15, 20, 30, 40, 50, 60, 80, 100, 150, 200, 300, 400, or more g/L. In some embodiments, the concentration of the carbohydrate is between about 50 and about 400 g/L, such as between about 100 and about 360 g/L, between about 120 and about 360 g/L, or between about 200 and about 300 g/L. In some embodiments, this concentration of carbohydrate includes the total amount of carbohydrate that is added before and/or during the culturing of the host cells.

[0260] Exemplary lipids are any substance containing one or more fatty acids that are C4 and above fatty acids that are saturated, unsaturated, or branched.

[0261] Exemplary fatty acids include compounds of the formula R-COOH, where "R" is a hydrocarbon. Exemplary unsaturated fatty acids include compounds where "R" includes at least one carbon-carbon double bond. Exemplary unsaturated fatty acids include, but are not limited to, oleic acid, vaccenic acid, linoleic acid, palmitelaidic acid, and arachidonic acid. Exemplary polyunsaturated fatty acids include compounds where "R" includes a plurality of carbon-carbon double bonds. Exemplary saturated fatty acids include compounds where "R" is a saturated aliphatic group. In some embodiments, the carbon source includes one or

more C₁₂-C₂₂ fatty acids, such as a C₁₂ saturated fatty acid, a C₁₄ saturated fatty acid, a C₁₆ saturated fatty acid, a C₁₈ saturated fatty acid, a C₂₀ saturated fatty acid, or a C₂₂ saturated fatty acid. In an exemplary embodiment, the fatty acid is palmitic acid. In some embodiments, the carbon source is a salt of a fatty acid (*e.g.*, an unsaturated fatty acid), a derivative of a fatty acid (*e.g.*, an unsaturated fatty acid), or a salt of a derivative of fatty acid (*e.g.*, an unsaturated fatty acid). Suitable salts include, but are not limited to, lithium salts, potassium salts, sodium salts, and the like. Di- and triglycerols are fatty acid esters of glycerol.

[0262] In some embodiments, the concentration of the lipid, fatty acid, monoglyceride, diglyceride, or triglyceride is at least or about 1 gram per liter of broth (g/L, wherein the volume of broth includes both the volume of the cell medium and the volume of the cells), such as at least or about 5, 10, 15, 20, 30, 40, 50, 60, 80, 100, 150, 200, 300, 400, or more g/L. In some embodiments, the concentration of the lipid, fatty acid, monoglyceride, diglyceride, or triglyceride is between about 10 and about 400 g/L, such as between about 25 and about 300 g/L, between about 60 and about 180 g/L, or between about 75 and about 150 g/L. In some embodiments, the concentration includes the total amount of the lipid, fatty acid, monoglyceride, diglyceride, or triglyceride that is added before and/or during the culturing of the host cells. In some embodiments, the carbon source includes both (i) a lipid, fatty acid, monoglyceride, diglyceride, or triglyceride and (ii) a carbohydrate, such as glucose. In some embodiments, the ratio of the lipid, fatty acid, monoglyceride, diglyceride, or triglyceride to the carbohydrate is about 1:1 on a carbon basis (*i.e.*, one carbon in the lipid, fatty acid, monoglyceride, diglyceride, or triglyceride per carbohydrate carbon). In particular embodiments, the amount of the lipid, fatty acid, monoglyceride, diglyceride, or triglyceride is between about 60 and 180 g/L, and the amount of the carbohydrate is between about 120 and 360 g/L.

[0263] Exemplary microbial polypeptide carbon sources include one or more polypeptides from yeast or bacteria. Exemplary plant polypeptide carbon sources include one or more polypeptides from soy, corn, canola, jatropha, palm, peanut, sunflower, coconut, mustard, rapeseed, cottonseed, palm kernel, olive, safflower, sesame, or linseed.

[0264] Exemplary renewable carbon sources include cheese whey permeate, cornsteep liquor, sugar beet molasses, barley malt, and components from any of the foregoing. Exemplary renewable carbon sources also include glucose, hexose, pentose and xylose present in biomass, such as corn, switchgrass, sugar cane, cell waste of fermentation processes, and protein by-product from the milling of soy, corn, or wheat. In some embodiments, the biomass carbon source is a lignocellulosic, hemicellulosic, or cellulosic

material such as, but are not limited to, a grass, wheat, wheat straw, bagasse, sugar cane bagasse, soft wood pulp, corn, corn cob or husk, corn kernel, fiber from corn kernels, corn stover, switch grass, rice hull product, or a by-product from wet or dry milling of grains (*e.g.*, corn, sorghum, rye, triticale, barley, wheat, and/or distillers grains). Exemplary cellulosic materials include wood, paper and pulp waste, herbaceous plants, and fruit pulp. In some embodiments, the carbon source includes any plant part, such as stems, grains, roots, or tubers. In some embodiments, all or part of any of the following plants are used as a carbon source: corn, wheat, rye, sorghum, triticale, rice, millet, barley, cassava, legumes, such as beans and peas, potatoes, sweet potatoes, bananas, sugarcane, and/or tapioca. In some embodiments, the carbon source is a biomass hydrolysate, such as a biomass hydrolysate that includes both xylose and glucose or that includes both sucrose and glucose.

[0265] In some embodiments, the renewable carbon source (such as biomass) is pretreated before it is added to the cell culture medium. In some embodiments, the pretreatment includes enzymatic pretreatment, chemical pretreatment, or a combination of both enzymatic and chemical pretreatment (*see*, for example, Farzaneh *et al.*, *Bioresource Technology* 96 (18): 2014-2018, 2005; U.S. Patent No. 6,176,176; U.S. Patent No. 6,106,888; which are each hereby incorporated by reference in their entireties, particularly with respect to the pretreatment of renewable carbon sources). In some embodiments, the renewable carbon source is partially or completely hydrolyzed before it is added to the cell culture medium.

[0266] In some embodiments, the renewable carbon source (such as corn stover) undergoes ammonia fiber expansion (AFEX) pretreatment before it is added to the cell culture medium (*see*, for example, Farzaneh *et al.*, *Bioresource Technology* 96 (18): 2014-2018, 2005). During AFEX pretreatment, a renewable carbon source is treated with liquid anhydrous ammonia at moderate temperatures (such as about 60 to about 100 °C) and high pressure (such as about 250 to about 300 psi) for about 5 minutes. Then, the pressure is rapidly released. In this process, the combined chemical and physical effects of lignin solubilization, hemicellulose hydrolysis, cellulose decrystallization, and increased surface area enables near complete enzymatic conversion of cellulose and hemicellulose to fermentable sugars. AFEX pretreatment has the advantage that nearly all of the ammonia can be recovered and reused, while the remaining serves as nitrogen source for microbes in downstream processes. Also, a wash stream is not required for AFEX pretreatment. Thus, dry matter recovery following the AFEX treatment is essentially 100%. AFEX is basically a dry-to-dry process. The treated renewable carbon source is stable for long periods and can be fed at very high solid loadings in enzymatic hydrolysis or fermentation processes. Cellulose and hemicellulose are well preserved in the AFEX process,

with little or no degradation. There is no need for neutralization prior to the enzymatic hydrolysis of a renewable carbon source that has undergone AFEX pretreatment. Enzymatic hydrolysis of AFEX-treated carbon sources produces clean sugar streams for subsequent fermentation use.

[0267] In some embodiments, the concentration of the carbon source (*e.g.*, a renewable carbon source) is equivalent to at least or about 0.1, 0.5, 1, 1.5, 2, 3, 4, 5, 10, 15, 20, 30, 40, or 50% glucose (w/v). The equivalent amount of glucose can be determined by using standard HPLC methods with glucose as a reference to measure the amount of glucose generated from the carbon source. In some embodiments, the concentration of the carbon source (*e.g.*, a renewable carbon source) is equivalent to between about 0.1 and about 20% glucose, such as between about 0.1 and about 10% glucose, between about 0.5 and about 10% glucose, between about 1 and about 10% glucose, between about 1 and about 5% glucose, or between about 1 and about 2% glucose.

[0268] In some embodiments, the carbon source includes yeast extract or one or more components of yeast extract. In some embodiments, the concentration of yeast extract is at least 1 gram of yeast extract per liter of broth (g/L, wherein the volume of broth includes both the volume of the cell medium and the volume of the cells), such as at least or about 5, 10, 15, 20, 30, 40, 50, 60, 80, 100, 150, 200, 300, or more g/L. In some embodiments, the concentration of yeast extract is between about 1 and about 300 g/L, such as between about 1 and about 200 g/L, between about 5 and about 200 g/L, between about 5 and about 100 g/L, or between about 5 and about 60 g/L. In some embodiments, the concentration includes the total amount of yeast extract that is added before and/or during the culturing of the host cells. In some embodiments, the carbon source includes both yeast extract (or one or more components thereof) and another carbon source, such as glucose. In some embodiments, the ratio of yeast extract to the other carbon source is about 1:5, about 1:10, or about 1:20 (w/w).

[0269] Additionally the carbon source may also be one-carbon substrates such as carbon dioxide, or methanol. Glycerol production from single carbon sources (*e.g.*, methanol, formaldehyde, or formate) has been reported in methylotrophic yeasts (Yamada *et al.*, *Agric. Biol. Chem.*, 53(2) 541-543, 1989, which is hereby incorporated by reference in its entirety, particularly with respect to carbon sources) and in bacteria (Hunter *et al.*, *Biochemistry*, 24, 4148-4155, 1985, which is hereby incorporated by reference in its entirety, particularly with respect to carbon sources). These organisms can assimilate single carbon compounds, ranging in oxidation state from methane to formate, and produce glycerol. The pathway of carbon assimilation can be through ribulose monophosphate, through serine, or through xylulose-

monophosphate (Gottschalk, Bacterial Metabolism, Second Edition, Springer-Verlag: New York, 1986, which is hereby incorporated by reference in its entirety, particularly with respect to carbon sources). The ribulose monophosphate pathway involves the condensation of formate with ribulose-5-phosphate to form a six carbon sugar that becomes fructose and eventually the three carbon product glyceraldehyde-3-phosphate. Likewise, the serine pathway assimilates the one-carbon compound into the glycolytic pathway via methylenetetrahydrofolate.

[0270] In addition to one and two carbon substrates, methylotrophic organisms are also known to utilize a number of other carbon containing compounds such as methylamine, glucosamine and a variety of amino acids for metabolic activity. For example, methylotrophic yeast are known to utilize the carbon from methylamine to form trehalose or glycerol (Bellion *et al.*, *Microb. Growth Cl Compd.*, [Int. Symp.], 7th ed., 415-32. Editors: Murrell *et al.*, Publisher: Intercept, Andover, UK, 1993, which is hereby incorporated by reference in its entirety, particularly with respect to carbon sources). Similarly, various species of *Candida* metabolize alanine or oleic acid (Sulter *et al.*, *Arch. Microbiol.* 153(5), 485-9, 1990, which is hereby incorporated by reference in its entirety, particularly with respect to carbon sources).

[0271] In some embodiments, cells are cultured in a standard medium containing physiological salts and nutrients (*see, e.g.*, Pourquie, J. *et al.*, *Biochemistry and Genetics of Cellulose Degradation*, eds. Aubert *et al.*, Academic Press, pp. 71-86, 1988; and Ilmen *et al.*, *Appl. Environ. Microbiol.* 63:1298-1306, 1997, hereby incorporated by reference, particularly with respect to cell media). Exemplary growth media are common commercially prepared media such as Luria Bertani (LB) broth, Sabouraud Dextrose (SD) broth, or Yeast medium (YM) broth. Other defined or synthetic growth media may also be used, and the appropriate medium for growth of particular host cells are known by someone skilled in the art of microbiology or fermentation science.

[0272] In addition to an appropriate carbon source, the cell medium desirably contains suitable minerals, salts, cofactors, buffers, and other components known to those skilled in the art suitable for the growth of the cultures or the enhancement of isoprene production (*see*, for example, WO 2004/033646 and references cited therein and WO 96/35796 and references cited therein, which are each hereby incorporated by reference in their entireties, particularly with respect cell medias and cell culture conditions). In some embodiments where an isoprene synthase, DXS, IDI, and/or MVA pathway nucleic acid is under the control of an inducible promoter, the inducing agent (*e.g.*, a sugar, metal salt or antimicrobial), is desirably added to the medium at a concentration effective to induce expression of an isoprene

synthase, DXS, IDI, and/or MVA pathway polypeptide. In some embodiments, cell medium has an antibiotic (such as kanamycin) that corresponds to the antibiotic resistance nucleic acid (such as a kanamycin resistance nucleic acid) on a vector that has one or more DXS, IDI, or MVA pathway nucleic acids.

Exemplary Production of Isoprene

[0273] In some embodiments, the cells are cultured in a culture medium under conditions permitting the production of isoprene by the cells. In some embodiments, the cells in culture produce isoprene at greater than or about 1, 10, 25, 50, 100, 150, 200, 250, 300, 400, 500, 600, 700, 800, 900, 1,000, 1,250, 1,500, 1,750, 2,000, 2,500, 3,000, 4,000, 5,000, or more nmole of isoprene/gram of cells for the wet weight of the cells/hour (nmole/g_{wcm}/hr). In some embodiments, the amount of isoprene is between about 2 to about 5,000 nmole/g_{wcm}/hr, such as between about 2 to about 100 nmole/g_{wcm}/hr, about 100 to about 500 nmole/g_{wcm}/hr, about 150 to about 500 nmole/g_{wcm}/hr, about 500 to about 1,000 nmole/g_{wcm}/hr, about 1,000 to about 2,000 nmole/g_{wcm}/hr, or about 2,000 to about 5,000 nmole/g_{wcm}/hr. The amount of isoprene in units of nmole/g_{wcm}/hr can be measured as disclosed in U.S. Patent No. 5,849,970, which is hereby incorporated by reference in its entirety, particularly with respect to the measurement of isoprene production. For example, two mL of headspace (*e.g.*, headspace from a culture such as 2 mL of culture cultured in sealed vials at 32 °C with shaking at 200 rpm for approximately 3 hours) are analyzed for isoprene using a standard gas chromatography system, such as a system operated isothermally (85 °C) with an n-octane/porasil C column (Alltech Associates, Inc., Deerfield, IL) and coupled to a RGD2 mercuric oxide reduction gas detector (Trace Analytical, Menlo Park, CA) (*see*, for example, Greenberg et al, *Atmos. Environ.* 27A: 2689-2692, 1993; Silver *et al.*, *Plant Physiol.* 97:1588-1591, 1991, which are each hereby incorporated by reference in their entireties, particularly with respect to the measurement of isoprene production). The gas chromatography area units are converted to nmol isoprene via a standard isoprene concentration calibration curve. In some embodiments, the value for the grams of cells for the wet weight of the cells is calculated by obtaining the A₆₀₀ value for a sample of the cell culture, and then converting the A₆₀₀ value to grams of cells based on a calibration curve of wet weights for cell cultures with a known A₆₀₀ value. In some embodiments, the grams of the cells is estimated by assuming that one liter of broth (including cell medium and cells) with an A₆₀₀ value of 1 has a wet cell weight of 1 gram. The value is also divided by the number of hours the culture has been incubating for, such as three hours.

[0274] In some embodiments, the cells in culture produce isoprene at greater than or about 1, 10, 25, 50, 100, 150, 200, 250, 300, 400, 500, 600, 700, 800, 900, 1,000, 1,250, 1,500, 1,750, 2,000, 2,500, 3,000, 4,000, 5,000, 10,000, 100,000, or more ng of isoprene/gram of cells for the wet weight of the cells/hr ($\text{ng/g}_{\text{wcm}}/\text{h}$). In some embodiments, the amount of isoprene is between about 2 to about 5,000 $\text{ng/g}_{\text{wcm}}/\text{h}$, such as between about 2 to about 100 $\text{ng/g}_{\text{wcm}}/\text{h}$, about 100 to about 500 $\text{ng/g}_{\text{wcm}}/\text{h}$, about 500 to about 1,000 $\text{ng/g}_{\text{wcm}}/\text{h}$, about 1,000 to about 2,000 $\text{ng/g}_{\text{wcm}}/\text{h}$, or about 2,000 to about 5,000 $\text{ng/g}_{\text{wcm}}/\text{h}$. The amount of isoprene in $\text{ng/g}_{\text{wcm}}/\text{h}$ can be calculated by multiplying the value for isoprene production in the units of $\text{nmole/g}_{\text{wcm}}/\text{hr}$ discussed above by 68.1 (as described in Equation 5 below).

[0275] In some embodiments, the cells in culture produce a cumulative titer (total amount) of isoprene at greater than or about 1, 10, 25, 50, 100, 150, 200, 250, 300, 400, 500, 600, 700, 800, 900, 1,000, 1,250, 1,500, 1,750, 2,000, 2,500, 3,000, 4,000, 5,000, 10,000, 50,000, 100,000, or more mg of isoprene/L of broth ($\text{mg/L}_{\text{broth}}$, wherein the volume of broth includes the volume of the cells and the cell medium). In some embodiments, the amount of isoprene is between about 2 to about 5,000 $\text{mg/L}_{\text{broth}}$, such as between about 2 to about 100 $\text{mg/L}_{\text{broth}}$, about 100 to about 500 $\text{mg/L}_{\text{broth}}$, about 500 to about 1,000 $\text{mg/L}_{\text{broth}}$, about 1,000 to about 2,000 $\text{mg/L}_{\text{broth}}$, or about 2,000 to about 5,000 $\text{mg/L}_{\text{broth}}$. The specific productivity of isoprene in mg of isoprene/L of headspace from shake flask or similar cultures can be measured by taking a 1 ml sample from the cell culture at an OD_{600} value of approximately 1.0, putting it in a 20 mL vial, incubating for 30 minutes, and then measuring the amount of isoprene in the headspace. If the OD_{600} value is not 1.0, then the measurement can be normalized to an OD_{600} value of 1.0 by dividing by the OD_{600} value. The value of mg isoprene/L headspace can be converted to $\text{mg/L}_{\text{broth}}/\text{hr}/\text{OD}_{600}$ of culture broth by multiplying by a factor of 38. The value in units of $\text{mg/L}_{\text{broth}}/\text{hr}/\text{OD}_{600}$ can be multiplied by the number of hours and the OD_{600} value to obtain the cumulative titer in units of mg of isoprene/L of broth.

[0276] The instantaneous isoprene production rate in $\text{mg/L}_{\text{broth}}/\text{hr}$ in a fermentor can be measured by taking a sample of the fermentor off-gas, analyzing it for the amount of isoprene (in units such as mg of isoprene per L_{gas}), and multiplying this value by the rate at which off-gas is passed through each liter of broth (*e.g.*, at 1 vvm (volume of air/volume of broth/minute) this is 60 L_{gas} per hour). Thus, an off-gas level of 1 mg/L_{gas} corresponds to an instantaneous production rate of 60 $\text{mg/L}_{\text{broth}}/\text{hr}$ at air flow of 1 vvm. If desired, the value in the units $\text{mg/L}_{\text{broth}}/\text{hr}$ can be divided by the OD_{600} value to obtain the specific rate in units of $\text{mg/L}_{\text{broth}}/\text{hr}/\text{OD}$. The average value of mg isoprene/ L_{gas} can be converted to the total product productivity (grams of isoprene per liter of fermentation broth, $\text{mg/L}_{\text{broth}}$) by multiplying this

average off-gas isoprene concentration by the total amount of off-gas sparged per liter of fermentation broth during the fermentation. Thus, an average off-gas isoprene concentration of 0.5 mg/L_{broth}/hr over 10 hours at 1 vvm corresponds to a total product concentration of 300 mg isoprene/L_{broth}.

[0277] In some embodiments, the cells in culture convert greater than or about 0.0015, 0.002, 0.005, 0.01, 0.02, 0.05, 0.1, 0.12, 0.14, 0.16, 0.2, 0.3, 0.4, 0.5, 0.6, 0.7, 0.8, 0.9, 1.0, 1.2, 1.4, or 1.6% of the carbon in the cell culture medium into isoprene. In some embodiments, the percent conversion of carbon into isoprene is between about 0.002 to about 1.6%, such as about 0.002 to about 0.005%, about 0.005 to about 0.01%, about 0.01 to about 0.05%, about 0.05 to about 0.15%, 0.15 to about 0.2%, about 0.2 to about 0.3%, about 0.3 to about 0.5%, about 0.5 to about 0.8%, about 0.8 to about 1.0%, or about 1.0 to about 1.6%. The percent conversion of carbon into isoprene (also referred to as “% carbon yield”) can be measured by dividing the moles carbon in the isoprene produced by the moles carbon in the carbon source (such as the moles of carbon in batched and fed glucose and yeast extract). This number is multiplied by 100% to give a percentage value (as indicated in Equation 1).

Equation 1

[0278] % Carbon Yield = (moles carbon in isoprene produced)/(moles carbon in carbon source) * 100

[0279] For this calculation, yeast extract can be assumed to contain 50% w/w carbon.

Equation 2

[0280] % Carbon Yield = (39.1 g isoprene * 1/68.1 mol/g * 5 C/mol)/[(181221 g glucose * 1/180 mol/g * 6 C/mol) + (17780 g yeast extract * 0.5 * 1/12 mol/g)] * 100 = 0.042%

[0281] One skilled in the art can readily convert the rates of isoprene production or amount of isoprene produced into any other units. Exemplary equations are listed below for interconverting between units.

Units for Rate of Isoprene production (total and specific)

Equation 3

[0282] 1 g isoprene/L_{broth}/hr = 14.7 mmol isoprene/L_{broth}/hr (total volumetric rate)

Equation 4

[0283] $1 \text{ nmol isoprene /g}_{\text{wcm}}/\text{hr} = 1 \text{ nmol isoprene /L}_{\text{broth}}/\text{hr}/\text{OD}_{600}$ (This conversion assumes that one liter of broth with an OD_{600} value of 1 has a wet cell weight of 1 gram.)

Equation 5

[0284] $1 \text{ nmol isoprene/g}_{\text{wcm}}/\text{hr} = 68.1 \text{ ng isoprene/g}_{\text{wcm}}/\text{hr}$ (given the molecular weight of isoprene)

Equation 6

[0285] $1 \text{ nmol isoprene/L}_{\text{gas}} \text{ O}_2/\text{hr} = 90 \text{ nmol isoprene/L}_{\text{broth}}/\text{hr}$ (at an O_2 flow rate of 90 L/hr per L of culture broth)

Equation 7

[0286] $1 \text{ }\mu\text{g isoprene/L}_{\text{gas}} \text{ isoprene in off-gas} = 60 \text{ }\mu\text{g isoprene/L}_{\text{broth}}/\text{hr}$ at a flow rate of 60 L_{gas} per L_{broth} (1 vvm)

Units for Titer (total and specific)

Equation 8

[0287] $1 \text{ nmol isoprene/mg cell protein} = 150 \text{ nmol isoprene/L}_{\text{broth}}/\text{OD}_{600}$ (This conversion assumes that one liter of broth with an OD_{600} value of 1 has a total cell protein of approximately 150 mg) (specific productivity)

Equation 9

[0288] $1 \text{ g isoprene/L}_{\text{broth}} = 14.7 \text{ mmol isoprene/L}_{\text{broth}}$ (total titer)

[0289] If desired, Equation 10 can be used to convert any of the units that include the wet weight of the cells into the corresponding units that include the dry weight of the cells.

Equation 10

[0290] Dry weight of cells = (wet weight of cells)/3.3

[0291] In some embodiments encompassed by the invention, a cell comprising a heterologous nucleic acid encoding an isoprene synthase polypeptide produces an amount of isoprene that is at least or about 2-fold, 3-fold, 5-fold, 10-fold, 25-fold, 50-fold, 100-fold, 150-fold, 200-fold, 400-fold, or greater than the amount of isoprene produced from a

corresponding cell grown under essentially the same conditions without the heterologous nucleic acid encoding the isoprene synthase polypeptide.

[0292] In some embodiments encompassed by the invention, a cell comprising a heterologous nucleic acid encoding an isoprene synthase polypeptide and one or more heterologous nucleic acids encoding a DXS, IDI, and/or MVA pathway polypeptide produces an amount of isoprene that is at least or about 2-fold, 3-fold, 5-fold, 10-fold, 25-fold, 50-fold, 100-fold, 150-fold, 200-fold, 400-fold, or greater than the amount of isoprene produced from a corresponding cell grown under essentially the same conditions without the heterologous nucleic acids.

Exemplary Isoprene Purification Methods

[0293] In some embodiments, any of the methods described herein further include recovering the isoprene. For example, the isoprene produced using the compositions and methods of the invention can be recovered using standard techniques, such as gas stripping, fractionation, adsorption/desorption, pervaporation, thermal or vacuum desorption of isoprene from a solid phase, or extraction of isoprene immobilized or absorbed to a solid phase with a solvent (see, for example, U.S. Patent Nos. 4,703,007 and 4,570,029, which are each hereby incorporated by reference in their entireties, particularly with respect to isoprene recovery and purification methods). In some embodiments, the recovery of isoprene involves the isolation of isoprene in a liquid form (such as a neat solution of isoprene or a solution of isoprene in a solvent). Gas stripping involves the removal of isoprene vapor from the fermentation off-gas stream in a continuous manner. Such removal can be achieved in several different ways including, but not limited to, adsorption to a solid phase, partition into a liquid phase, or direct condensation. In some embodiments, membrane enrichment of a dilute isoprene vapor stream above the dew point of the vapor resulting in the condensation of liquid isoprene.

[0294] The recovery of isoprene may involve one step or multiple steps. In some embodiments, the removal of isoprene vapor from the fermentation off-gas and the conversion of isoprene to a liquid phase are performed simultaneously. For example, isoprene can be directly condensed from the off-gas stream to form a liquid. In some embodiments, the removal of isoprene vapor from the fermentation off-gas and the conversion of isoprene to a liquid phase are performed sequentially. For example, isoprene may be adsorbed to a solid phase and then extracted from the solid phase with a solvent.

[0295] In some embodiments, any of the methods described herein further include purifying the isoprene. For example, the isoprene produced using the compositions and methods of the invention can be purified using standard techniques. Purification refers to a process through which isoprene is separated from one or more components that are present when the isoprene is produced. In some embodiments, the isoprene is obtained as a substantially pure liquid. Examples of purification methods include (i) distillation from a solution in a liquid extractant and (ii) chromatography. As used herein, “purified isoprene” means isoprene that has been separated from one or more components that are present when the isoprene is produced. In some embodiments, the isoprene is at least about 20%, by weight, free from other components that are present when the isoprene is produced. In various embodiments, the isoprene is at least or about 25%, 30%, 40%, 50%, 60%, 70%, 75%, 80%, 90%, 95%, or 99%, by weight, pure. Purity can be assayed by any appropriate method, *e.g.*, by column chromatography, HPLC analysis, or GC-MS analysis.

[0296]

[0297] Crystal Structure of Isoprene Synthase

[0298] The invention also contemplates crystalline forms of plant isoprene synthase (*e.g.*, poplar and kudzu) and its variants as described *supra* and in the Examples. In one embodiment, the invention comprises any polypeptide which has the crystal structure of poplar isoprene synthase as disclosed in Table 16-7.

EXPERIMENTAL

[0299] The following examples are provided in order to demonstrate and further illustrate certain preferred embodiments and aspects of the present invention and are not to be construed as limiting the scope thereof.

[0300] In the experimental disclosure which follows, the following abbreviations apply: °C (degrees Centigrade); rpm (revolutions per minute); H₂O (water); diH₂O (deionized water); aa and AA (amino acid); bp (base pair); kb (kilobase pair); kD (kilodaltons); gm (grams); µg and ug (micrograms); mg (milligrams); ng (nanograms); µl and ul (microliters); ml (milliliters); mm (millimeters); qs (quantity sufficient); nm (nanometers); µm and um (micrometer); M (molar); mM (millimolar); µM and uM (micromolar); pM (picomolar); U (units); MW (molecular weight); sec (seconds); min (minute/minutes); hr (hour/hours); OD₆₀₀ (optical density at 600 nm); BSA (bovine serum albumin); DMAPP (dimethylallyl diphosphate); DTT (dithiothreitol); EtOH (ethanol); IPTG (isopropyl-beta-D-

thiogalactopyranoside); isoprene (2-methyl-1,3-butadiene); IspS (isoprene synthase); PAGE (polyacrylamide gel electrophoresis); PBS (phosphate buffered saline [150 mM NaCl, 10 mM sodium phosphate buffer, pH 7.2]); and SDS (sodium dodecyl sulfate).

[0301] The following abbreviations apply to companies whose products or services may have been referred to in the experimental examples: Agilent (Agilent Technologies, Santa Clara, CA); Becton Coulter (Becton Coulter, Inc., Fullerton, CA); Bio-Rad (Bio-Rad Laboratories, Hercules, CA); Cayman Chemical (Cayman Chemical Co., Ann Arbor, MI); CTC Analytics (CTC Analytics A.G., Zwingen, Switzerland); EMS (Electron Microscopy Supply, Hatfield, PA); Epicentre (Epicentre Biotechnologies, Madison, WI); Integrated DNA Technologies (Integrated DNA Technologies, Coralville, IA); Invitrogen (Invitrogen Corp., Carlsbad, CA); Molecular Dynamics (Molecular Dynamics, Sunnyvale, CA); Novagen (Novagen, Inc., Madison, WI); Perkin Elmer (Perkin Elmer, Waltham, MA); Roche (Roche Applied Science, Indianapolis, IN); Sigma (Sigma-Aldrich, St. Louis, MO); Stratagene (Stratagene Cloning Systems, La Jolla, CA); Qiagen (Qiagen, Inc., Valencia, CA); Takara (Takara Bio USA, Madison, WI); Thomson Instrument (Thomson Instrument Co., Oceanside, CA); V&P Scientific (V&P Scientific, Inc., San Diego, CA); and Zinsser (Zinsser North America, Northridge, CA).

EXAMPLE 1

Cloning of Kudzu Isoprene Synthase for Expression in Recombinant Bacteria

[0302] In this example, methods used to produce kudzu isoprene synthase (IspS) in *E. coli* are described. The protein sequence for the kudzu (*Pueraria montana*) isoprene synthase gene (IspS) was obtained from GenBank (AAQ84170). A kudzu isoprene synthase gene, optimized for *E. coli* codon usage, was purchased from DNA2.0 (Menlo Park, CA), and is set forth as SEQ ID NO:1 (Figure 1). The isoprene synthase gene was removed from the supplied plasmid by restriction endonuclease digestion with *Bsp*LU11I/*Pst*I, gel-purified, and ligated into pTrcHis2B (Invitrogen) that had been digested with *Nco*I/*Pst*I. The construct was designed such that the stop codon in the isoprene synthase gene was 5' to the *Pst*I site. As a result, when the construct was expressed the His-Tag is not attached to the isoprene synthase protein. The resulting plasmid, pTrcKudzu, was verified by sequencing.

[0303] The isoprene synthase gene was also cloned into pET16b (Novagen). In this case, the isoprene synthase gene was inserted into pET16b such that the recombinant isoprene synthase protein contained the N-terminal His tag. The isoprene synthase gene was amplified

from pTrcKudzu by PCR using the primer set pET-His-Kudzu-2F: 5'-CGTGAGATCA TATGTGTGCG ACCTCTTCTC AATTTAC (SEQ ID NO:3) and pET-His-Kudzu-R: 5'-CGGTCGACGG ATCCCTGCAG TTAGACATAC ATCAGCTG (SEQ ID NO:4). These primers added an *NdeI* site at the 5'-end and a *BamHI* site at the 3' end of the gene respectively. The plasmid pTrcKudzu, described above, was used as template DNA, HERCULASE DNA polymerase (Stratagene) was used according to manufacturer's directions, and primers were added at a concentration of 10 pM. The PCR was carried out in a total volume of 25 µl. The PCR product was digested with *NdeI/BamHI* and cloned into pET16b digested with the same enzymes. The ligation mix was transformed into *E. coli* Top10 (Invitrogen) and the correct clone selected by sequencing. The resulting plasmid designated pETNHiskudzu is then transformed into BL21(λDE3)pLysS (Novagen) cells for expression from the T7 promoter.

[0304] The kudzu isoprene synthase gene was also cloned into the low copy number plasmid pCL1920 (Lerner and Inouye, Nucl Acids Res, 18:4631, 1990). Primers were used to amplify the kudzu isoprene synthase gene from pTrcKudzu described above. The forward primer added a *HindIII* site and an *E. coli* consensus RBS to the 5' end. The *PstI* cloning site was already present in pTrcKudzu just 3' of the stop codon so the reverse primer was constructed such that the final PCR product includes the *PstI* site. The sequences of the primers were: *HindIII*-rbs-Kudzu F: 5'-CATATGAAAG CTTGTATCGA TTAAATAAGG AGGAATAAAC C (SEQ ID NO:5) and *BamHI*-Kudzu R: 5'-CGGTCGACGG ATCCCTGCAG TTAGACATAC ATCAGCTG (SEQ ID NO:4). The PCR product was amplified using HERCULASE DNA polymerase (Stratagene) with primers at a concentration of 10 pM and with 1 ng of template DNA (pTrcKudzu). The amplification protocol included 30 cycles of (95° C for 1 minute, 60° C for 1 minute, 72° C for 2 minutes). The product was digested with *HindIII* and *PstI* and ligated into pCL1920, which had also been digested with *HindIII* and *PstI*. The ligation mix was transformed into *E. coli* Top10. Several transformants were verified by sequence analysis. The resulting plasmid was designated pCL-lac-Kudzu.

[0305] In order to remove the beta-lactamase gene, pTrcKudzu was digested with *BspHI*, treated with shrimp alkaline phosphatase (SAP), incubated at 65°C for 10 min to heat kill the SAP, then end-filled by incubating with 2 units of Klenow fragment (New England BioLabs) and dNTPS. The 5 kb fragment was purified from an agarose gel and ligated to the Kan(R) gene. The Kan(R) gene was prepared by PCR amplification from pCR-Blunt-II-TOPO (Invitrogen) using primers MCM22 and MCM23 and Taq DNA polymerase according

to the Manufacturer's instructions. The PCR fragment was digested with *Hind*III and *Pvu*I, and end-filled using Klenow Fragment and dNTPs. The ligation mixture was transformed into *E. coli* Top 10 chemically competent cells and a transformant carrying a plasmid conferring kanamycin resistance, pTrcKudzu(kan), was selected on Luria Agar containing kanamycin (50µg/ml). The sequences of the primers were: MCM22 5'- gatcaagctt AACCGGAATTGCCAGCTG (SEQ ID NO:15); and MCM23 5'- gatccgatcgTCAGAAGAACTCGTCAAGAAGGC (SEQ ID NO:16).

EXAMPLE 2

Cloning of Poplar Isoprene Synthase for Expression in Recombinant Bacteria

[0306] In this example, methods used to produce poplar isoprene synthase (IspS) in *E. coli* are described. The protein sequence for the poplar (*Populus alba* x *Populus tremula*) isoprene synthase (Schnitzler et al., Planta 222:777-786, 2005) was obtained from GenBank (CAC35696). A gene, codon optimized for *E. coli*, was purchased from DNA2.0 and is set forth as SEQ ID NO:6 (Figure 3). The isoprene synthase gene was removed from the supplied plasmid (p9796-poplar) by restriction endonuclease digestion with *Bsp*LU11I/*Pst*I, gel-purified, and ligated into pTrcHis2B that had been digested with *Nco*I/*Pst*I. The construct is cloned such that the stop codon in the insert is before the *Pst*I site, which results in a construct in which the His-Tag is not attached to the isoprene synthase protein. The resulting plasmid, pTrcPoplar, was verified by sequencing using commercially available primers that hybridize within the vector sequence (Forward and Reverse), as well as the primer Poplar InSeq 5'-GAGAAAATCG GTAAGGAACT GG (SEQ ID NO:8).

EXAMPLE 3

Isoprene Production in Recombinant Bacteria

[0307] In this example, methods used to produce and measure isoprene in recombinant *E. coli* are described.

I. Determination of isoprene production

[0308] For the shake flask cultures, one ml of a culture was transferred from shake flasks to 20 ml CTC headspace vials (Agilent vial Catalog No. 5188 2753, and cap Catalog No. 5188 2759). The cap was screwed on tightly and the vials incubated at the equivalent temperature

with shaking at 250 rpm. After 30 minutes, the vials were removed from the incubator and analyzed as described below. In cases where isoprene production in fermentors was determined, samples were taken from the off-gas of the fermentor and analyzed directly.

[0309] The analysis was performed using an Agilent 6890 GC/MS system interfaced with a CTC Analytics CombiPAL autosampler operating in headspace mode. An Agilent HP-5MS GC/MS column (30 m x 0.25 mm; 0.25 μ m film thickness) was used for separation of analytes. The sampler was set up to inject 500 μ L of headspace gas. The GC/MS method utilized helium as the carrier gas at a flow of 1 ml/min. The injection port was held at 250° C with a split ratio of 50:1. The oven temperature was held at 37° C for the 2 min duration of the analysis. The Agilent 5793N mass selective detector was run in single ion monitoring (SIM) mode on m/z 67. The detector was switched off from 1.4 to 1.7 minutes to allow the elution of permanent gases. Under these conditions, isoprene (2-methyl-1,3-butadiene) was observed to elute at 1.78 minutes. A calibration table was used to quantify the absolute amount of isoprene and was found to be linear from 1 μ g/L to 200 μ g/L. The limit of detection was estimated to be 50 to 100 ng/L using this method.

II. Production of isoprene in shake flasks

[0310] The vectors described above were introduced into *E. coli* strain BL21(λ DE3)pLysS (Novagen) to produce strains BL21/ptrcKudzu, BL21/pCL-lac-Kudzu and BL21/pETHisKudzu. The strains were spread for isolation onto LA (Luria agar) containing the appropriate antibiotic (50 μ g/ml carbenicillin for BL21/ptrcKudzu and BL21/pETHisKudzu or 50 μ g/ml spectinomycin for BL21/pCL-lac-Kudzu) and incubated overnight at 37° C. Single colonies were inoculated into 250 ml baffled shake flasks containing 20 ml Luria Bertani broth (LB) and the appropriate antibiotic. Cultures were grown overnight at 20° C with shaking at 200 rpm. The OD₆₀₀ of the overnight cultures was measured and the cultures were diluted into a 250 ml baffled shake flask containing 30 ml MAGICMEDIA expression medium (Invitrogen) containing the appropriate antibiotic to an OD₆₀₀ ~ 0.05. The culture was incubated at 30° C with shaking at 200 rpm. When the OD₆₀₀ ~ 0.5 - 0.8, 400 μ M IPTG was added and the cells were incubated for a further 6 hours at 30° C with shaking at 200 rpm. At 0, 2, 4 and 6 hours after induction with IPTG, 1 ml aliquots of the cultures were collected, the OD₆₀₀ was determined and the amount of isoprene produced was measured as described above.

III. Production of Isoprene from BL21/ptrcKudzu in 14 liter fermentation

[0311] Large-scale production of isoprene from *E. coli* containing the recombinant kudzu isoprene synthase gene was determined from a fed-batch culture. The recipe for the fermentation media (TM2) per liter of fermentation medium was as follows: K_2HPO_4 13.6 g, KH_2PO_4 13.6 g, $MgSO_4 \cdot 7H_2O$ 2 g, citric acid monohydrate 2 g, ferric ammonium citrate 0.3 g, $(NH_4)_2SO_4$ 3.2 g, yeast extract 5 g, 1000X Modified Trace Metal Solution 1 ml. All of the components were added together and dissolved in diH_2O . The pH was adjusted to 6.8 with potassium hydroxide (KOH) and qs to volume. The final product was filter sterilized with 0.22 μ M filter, but not autoclaved. The recipe for 1000X Modified Trace Metal Solution was as follows: Citric Acids * H_2O 40 g, $MnSO_4 \cdot H_2O$ 30 g, NaCl 10 g, $FeSO_4 \cdot 7H_2O$ 1 g, $CoCl_2 \cdot 6H_2O$ 1 g, $ZnSO_4 \cdot 7H_2O$ 1 g, $CuSO_4 \cdot 5H_2O$ 100 mg, H_3BO_3 100 mg, $NaMoO_4 \cdot 2H_2O$ 100 mg. Each component was dissolved one at a time in diH_2O , the pH was adjusted to 3.0 with HCl/NaOH, then qs to volume and filter sterilized with a 0.22 μ filter.

[0312] This experiment was carried out in 14 L bioreactor to monitor isoprene formation from glucose at the desired fermentation, pH 6.7 and temperature 34° C. An inoculum of *E. coli* strain BL21/ptrcKudzu taken from a frozen vial was prepared in soytone-yeast extract-glucose medium in two 600 ml flasks.. After the inoculum grew to $OD_{550} = 0.6$, two 600 ml flasks were centrifuged and the contents resuspended in 70 ml supernatant to transfer the cell pellet (70 ml of OD 3.1 material) to the bioreactor. At various times after inoculation, samples were removed and the amount of isoprene produced was determined as described above.

EXAMPLE 4

Selection of Sites for Improvement of Plant Isoprene Synthases

[0313] The isoprene synthases of plants were expected to be homologous to the terpene synthases. The three-dimensional structures of two homologous terpene synthases have been determined from bornyl diphosphate synthase (pdb entry 1N1B) and 5-epi-aristolochene synthase (pdb entry 5EAU). These enzymes share only 32% homology but their tertiary structure is conserved. In addition, the structures of intermediate complexes with both enzymes have shown that not only tertiary folding, but also detailed interactions in the active sites of these enzymes are highly conserved.

[0314] The kudzu and poplar isoprene synthases have higher sequence identity than was seen between the bornyl diphosphate synthase and the 5-epi-aristolochene synthase as shown in Table 4-1 below.

Table 4-1. Percent Identity of Various Enzymes

	BDP-synthase	5EA-synthase	Kudzu IspS
Poplar IspS	40.1	32.9	54.4
Kudzu IspS	40.7	33.8	
5ES synthase	31.9		

[0315] A homology model of the poplar isoprene synthase has been made based on the bornyl diphosphate synthase (BDP-synthase) pdb entry 1N24 (~40% sequence identity). The homology model appears to be plausible based on the close similarity of 10 trial models created using the program MOE written and supported by The Chemical Computing Group, Inc. The plausibility is based on the conservation of common amino acid residues at sites found to be involved in catalysis in the BDP-synthase structure.

[0316] A comparison of the active site from the structure of BDP-synthase and the homology model of poplar IspS indicates that the active site involved in metal ion binding and phosphate recognition is conserved. In particular, Lys 272, Asp 309, Asp 313, Glu 387, Arg 450 and Asn 453 of poplar IspS were observed to overlap equivalent residues in BDP-synthase. In this example, amino acid residue positions for poplar IspS are derived from SEQ ID NO:7. The positioning of an intermediate of the BDP-synthase was also compared with the poplar IspS homology model. Based on this, it was possible to identify the analogous binding region and the approach direction that isopentenyl diphosphate would require in order to bind and react with the poplar IspS enzyme.

[0317] A homology model of the kudzu isoprene synthase has been made based on bornyl diphosphate synthase pdb entry 1N24 having (~40% sequence identity). A comparison of the active site from the structure of BDP-synthase with the homology model of kudzu IspS indicates that numerous active site residues involved in metal ion binding and phosphate recognition are conserved. In particular, Arg 269, Asp 306, Asp 310, Glu 384, Arg 450 and Asn 453 of kudzu IspS were observed to overlap equivalent residues in BDP-synthase.

[0318] A comparison of the active site residues identified in the homology models of poplar and kudzu IspS revealed that residues from one homology model are also quite homologous with similar residues, appearing with only minor shifts in the relative position numbers for some of the residues, in the other homology model. Based on the homology models, sites in poplar and kudzu IspS were identified as candidates for mutagenesis to produce variant IspS enzymes with improved performance. Briefly, sites were selected in the

IspS that might alter the interaction of the metal binding, the diphosphate recognition, the IPP chain binding and/or the approach to the active site.

I. Diphosphate/ Metal Binding Sites

[0319] The side chains of amino acid residues in the poplar IspS that are found in proximity to the metal and diphosphate (DPP) binding side chains were identified. These residues include Phe 384, Tyr 402, Ala 406, Ser 409, Ala 460 and Asn 469. The inventors note that Lys 272 is incorrect based on homology to other known poplar IspS sequences, which have an Arg at this position.

II. Substrate Access Loop

[0320] The substrate access loop of poplar IspS is in a region that deviates from the BDP-synthase structure. In the BDP-synthase structure the residues form a segment that creates a cover. Without being bound by theory, the inventors expect that this segment in the actual three-dimensional structure of poplar IspS will form a similar structure. As such the residues in this loop, including residues 455-466, will be in a position to alter the activity of the poplar IspS enzyme. In the poplar IspS enzyme residues 454-466 have the following sequence: LASASAEIARGET (SEQ ID NO:9).

III. Isoprenyl Binding Site

[0321] The complex of BDP-synthase and the product of the reaction, bornyl diphosphate (pdb entry 1N24), was used to identify residues in the poplar model that may modulate substrate specificity and/or reaction rate (altered on and off rates of substrate and product). These residues include Arg 274, Trp 281 Phe 302, Val 305, Ser 411, Gln 415, Phe 449, Ser 537 and Glu 540.

Table 4-2. Candidate Mutagenesis Sites

	Poplar Amino Acid	Kudzu Amino Acid / codon
DPP/metal sites	Phe 384	Phe 381 / 1141-1143
	Tyr 402	Tyr 399 / 1195-1197
	Ala 406	Ala 403 / 1207-1209
	Ser 409	Ser 406 / 1216-1218
	Asn 469	Asn 469 / 1405-1407

Extra DPP sites	Tyr 312	Tyr 309 / 925-927
	Asp 313	Asp 310 / 928-930
	Leu 380	Leu 377 / 1129-1131
	Glu 387	Glu 384 / 1150-1152
	Asn 404	Asn 402 / 1204-1206
	Ser 410	Ser 407 / 1219-1221
N_terminal access segment	22	20 / 58-60
	23	21 / 61-63
	24	22 / 64-66
	25	23 / 67-69
	26	24 / 70-72
	27	25 / 73-75
Substrate access loop	Leu 454	Ala 456 / 4102-4104
	Ala 455	Thr 457 / 4105-4107
	Ser 456	Ser 458 / 4108-4110
	Ala 457	Ala 459 / 4111-4113
	Ser 458	Ala 460 / 4114-4116
	Ala 459	Glu 461 / 4117-4119
	Glu 460	Leu 462 / 4120-4122
	Ile 461	Glu 463 / 4123-4125
	Ala 462	
	Arg 463	Arg 464 / 4126-4128
	Gly 464	Gly 465 / 4129-4131
	Glu 465	Glu 466 / 4132-4134
	Thr 466	Thr 467 / 4135-4137
		Thr 468 / 4138-4140
Isoprenyl binding site	Arg 274	Arg 271 / 811-813
	Trp 281	Trp 278 / 832-834
	Phe 302	Phe 299 / 895-897
	Val 305	Val 302 / 904-906
	Ser 411	Ser 408 / 1222-1224
	Gln 415	
	Phe 449	Phe 449 / 1345-1347
	Ser 537	Ser 458 / 1372-1374
	Glu 540	Tyr 531 / 1591-1593
Extra Isoprenyl sites	Gly 412	Gly 409 / 1225-1227

	Leu 414	Ala 411 / 1231-1233
	Leu 416	Leu 413 / 1237-1239
	Leu 521	Met 523 / 1567-1569
	Ser 525	Ser 527 / 1579-1581

EXAMPLE 5

Mutation of Non-Conserved Cysteines in Kudzu Isoprene Synthase

[0322] The kudzu and poplar isoprene synthase (IspS) homology models, based on the bornyl diphosphate synthase crystal structure, were compared with respect to the positions of the cysteine residues. Cysteines have the potential to form disulfide bonds and stabilize structures. The non-conserved cysteines, contemplated to affect solubility and/or activity, were altered by site-directed mutagenesis. The kudzu IspS amino acid sequence used for the modeling is shown in Figure 2 (SEQ ID NO:2). There are eight cysteines in kudzu IspS at positions 57, 291, 343, 378, 421, 446, 451 and 529 (relative to the mature form of the protein) as shown in the homology model of Figure 5. In contrast, there are five cysteines in poplar IspS amino acid sequence has five cysteines, as shown in the homology model of Figure 6 upon which the kudzu cysteines are superimposed. Several of the cysteines are apparently conserved between the poplar and kudzu IspS sequences indicating that these positions are important in stabilizing the structure, activity and/or other protein function. The remaining cysteines (nonconserved residues 57, 291, 421 and 446) in kudzu were mutated to serine as described herein.

I. Mutagenesis

[0323] The QUIKCHANGE® Multi-Site Directed Mutagenesis Kit (Stratagene) was used as per the manufacturer's directions. The following primers were utilized for mutagenesis:

C57S-F 5'- CTGGAGGAAGAAGTTCGC TCC
 ATGATCAACCGTGTAGAC (SEQ ID NO:10);
 C291S-F 5'- CGCCAGACCCGCAGTTTGGTGAA TCT
 CGCAAAGCTGTTACTAAAATG (SEQ ID NO:11);
 C421S-F 5'-CGCCGTCTTACTTTTCCGTA TCC
 CAGCAGCAGGAAGACATC (SEQ ID NO:12);
 C446S-F 5'-CATGGTCTGGTGCGTTCTAGC TCC
 GTTATCTTCCGCCTGTGC (SEQ ID NO:13); and
 C421S-R 5'-GATGTCTTCCTGCTGCTG GGA

TACGGAAAAGTAAGACGGCG (SEQ ID NO:14).

[0324] The plasmid pTrcKudzu(kan) described in Example 1 was used as template DNA. The primers C57S-F, C291S-F, C421S-F, and C446-F were combined in a single reaction (100 pmol). Template DNA was added (~200 nanograms) and 0.5 μ l of Quiksolution was added to the recommended volumes of enzyme and buffer. The PCR reaction was carried out in an Eppendorf PCR machine using an annealing temperature of 55°C and an extension time of 12 minutes for 30 cycles. Other parameters of the cycle were as indicated in the instructions. The PCR mix was treated with *DpnI* for 4 hours at 37°C (2 x 1 μ l for 2 h each) and then 5 μ l of the reaction were transformed into *E. coli* Top10 (Invitrogen) chemically competent cells and plated on Luria agar containing kanamycin (50 μ g/ml). After overnight incubation at 37°C, several colonies were picked and inoculated into 5 ml of Luria Broth containing kanamycin (50 μ g/ml). The plasmids were isolated using the QIAprep Spin Miniprep kit (Qiagen), and the IspS genes were sequenced in their entirety. Various single and combinations of mutations were made as indicated in the Table 5-1 below.

Table 5-1. BL21(λ DE3) Cells Transformed with Mutated pTrcKudzu Plasmids*

Strain*	C57S	C291S	C421S	C446S
C1	+	+		
C2		+		+
C4	+	+		+
C6	+			+
C11				+
C20			+	
C6-4	+		+	+

[0325] All the variant plasmids were transformed into chemically competent BL21(λ DE3) cells (Novagen). In a second reaction pTrcKudzu(kan) and plasmid DNA isolated from C6 were used as templates in a QUIKCHANGE® site directed mutagenesis kit (Stratagene) single site reaction with C421S-F and C421S-R primers. After confirmation by sequencing, two additional strains were obtained.

II. Cell Growth and Isoprene Production

[0326] Cells were grown in 5 ml tubes containing Luria Broth supplemented with 50 mg/L kanamycin at 30°C overnight with agitation. These cultures were diluted into TM3 broth supplemented with 10 g/L glucose and 50 mg/L kanamycin. The culture volume was 25 ml in a 250 ml baffled Bellco Delong flask in which cells were grown at 30°C with agitation (225 rpm). Samples were taken aseptically, as indicated, for optical density measurements at A₆₀₀. The results are shown in Figure 7. The cultures were induced at 3.33 hrs with 200 µM IPTG and allowed to continue growth until harvest at 7.8 hr. The cultures were centrifuged at 10,000 X g for 10 min, the supernatants decanted and the cell pellets frozen at -80°C overnight.

[0327] Frozen cell pellets were thawed and resuspended in 2 ml PEB (50 mM Tris-HCl, pH 8.0, 20 mM MgCl₂, 2 mM dithiothreitol, and 50% [v/v] glycerol). Cells were lysed by French pressure cell disruption, one pass, at 20,000 psi. The lysate was centrifuged for 15 min at 10,000 X g. The supernatants were decanted and the pellets resuspended in 2 ml PEB. The pellets and supernatants were analyzed by SDS-PAGE, 4-12% NuPage gels (Invitrogen), run in MES buffer under reducing conditions. The molecular weight standard was SeeBlue2 (Invitrogen). The results are shown in Figure 8. The IspS protein concentrations were estimated using the BCA assay (Pierce) using BSA as standard (Table 5-2).

III. Assays for Isoprene Synthase Activity and Solubility

[0328] Briefly the activity of the supernatants was measured by reaction with DMAPP, and the isoprene evolved was quantified by GC/MS.

[0329] Headspace Assay. A sample of 200 µl of the desired culture is inoculated into 2 ml CTC headspace vials (Agilent vial Catalog No. 5188 2753, and cap Catalog No. 5188 2759). The cap was screwed on tightly and the vials were incubated at 37°C with shaking at 250 rpm. After 30 minutes the vials were removed from the incubator and cooled briefly with ambient tap water. The vials were placed into the CombiPal Headspace auto sampler for analysis by GC-MS. The analysis was performed using an Agilent 6890 GC/MS system interfaced with a CTC Analytics CombiPAL autosampler operating in headspace mode. An Agilent HP-5MS GC/MS column (30 m x 0.25 mm; 0.25 µm film thickness) was used for separation of analytes. The sampler was set up to inject 500 µL of headspace gas. The GC/MS method utilized helium as the carrier gas at a flow of 1 ml/min. The injection port was held at 250°C with a split ratio of 50:1. The oven temperature was held at 37°C for the 2 min duration of the analysis. The Agilent 5793N mass selective detector was run in single ion

monitoring (SIM) mode on m/z 67. The detector was switched off from 1.4 to 1.7 minutes to allow the elution of permanent gases. Under these conditions isoprene (2-methyl-1,3-butadiene) was observed to elute at 1.78 minutes. A calibration table was used to quantify the absolute amount of isoprene and was found to be linear from 1 $\mu\text{g/L}$ to 200 $\mu\text{g/L}$. The limit of detection was estimated to be 50-100 ng/L using this method.

[0330] DMAPP Assay. An aliquot of 95 μl of the supernatant fraction from the centrifuged French Pressure cell lysate was added to the headspace vials. A 5 μL aliquot of 8 mM DMAPP in 100 mM potassium phosphate buffer, pH 8.2 was added, the vials sealed and allowed to incubate at room temperature for 30 min. The amount of isoprene produced was measured by GC/MS as described above and reported in Table 5-2.

Table 5-2. Isoprene Synthase Activity from Crude Extracts of Cysteine Mutants

Supernatant	Activity (mU/ml)	Protein (mg/ml)	Specific Activity (mU/mgP)
Wt	59.0	2.48	23.79
C1	5.3	2.64	2.01
C2	0.2	3.52	0.06
C4	0.5	3.22	0.16
C6	9.0	3.32	2.71
C11	0.1	4.26	0.03

[0331] The values shown in Table 5-2 are averages of reactions with two different concentrations of extract. All proteins containing any of the cysteine mutations resulted in severe diminution of enzyme activity and an apparent decrease of soluble protein as judged by a relative increase in the proportion of protein in the insoluble (pellet) fraction.

EXAMPLE 6

Mutation of Residues in Poplar Isoprene Synthase

[0332] Alignment of the amino acid sequences of kudzu and poplar isoprene synthases with other synthases was done using Vector Nti (Invitrogen). The aligned sequences included: beta-ocimene synthase *Lotus corniculatus* (AAT86042); putative terpene synthetase *Medicago trunculata* (AAV36465); hypothetical protein *Vitis vinifera* (can65805); hypothetical protein *Vitis vinifera* (CAN62729); pinene synthase *Quercus ilex* (CAK55186); IspS *Pueraria montana* (kudzu) Sharkey (AAQ84170); monoterpene synthase *Eucalyptus globulus* (BAF02832); IspS *Populus nigra* Fortunati (CAL69918); IspS *Populus tremuloides*

Sharkey (AAQ16588); IspS *Populus alba* (BAD98243); and IspS *Populus alba x tremula* Zimmer (CAC35696). The sequence from the database of *Populus alba x tremula* (CAC35696) exhibited different amino acids at positions 272 and 497 that were otherwise highly conserved. Additionally based on analysis of the homology model of poplar IspS, position 453 was identified as a third candidate for mutagenesis.

I. Mutagenesis

[0333] The QUIKCHANGE® Multi-Site Directed Mutagenesis kit (Stratagene) was used as per the manufacturer's directions to introduce the following mutations singly and in combination into the *Populus alba x tremula* IspS sequence (SEQ ID NO:7): K272R; C497W; and N453D. The following primers were utilized for mutagenesis:

Poplar K272R 5'-ccaaactgcacttcg ctcgtgaccgcctgattgag (SEQ ID NO:17);

Poplar N453D 5'-atctttcgctgtgacgacacctggcaagc (SEQ ID NO:18); and

Poplar C497W 5'-tgaatctgatcgacgaaacctggaagaaaatgaacaaagaaaaac (SEQ ID NO:19).

The following primer, Poplar InSeq, 5'-gagaaaatcggttaaggaactgg (SEQ ID NO:8) was used for sequencing.

[0334] Mutagenesis was done according to the manufacturer's directions, with all three mutagenesis primers being added to a single reaction mix (100 ng each) with pTrcPoplar as the template DNA (100 ng). Addition of 0.5 µl of Quik Solution aided the mutagenesis reaction. The suggested PCR cycle was run with an annealing temperature of 55°C and an extension time of 12 min. Other parameters were as indicated in the instructions. The PCR mix was digested with *DpnI* for 4 hrs at 37°C (1 µl each x 2 h) and then 5 µl of the reaction was transformed into *E. coli* Top10 chemically competent cells (Invitrogen). Several colonies were selected and grown up in 5 ml of Luria Broth (LB) containing carbenicillin (50 µg/ml). Plasmids were isolated using the Qiagen QIAprep spin mini-prep kit and sent for sequencing using forward and reverse primers that hybridized to the vector, as well as the Poplar InSeq primer.

[0335] The resulting variants of pTrcPoplar were obtained: pTrcPoplar K272R, pTrcPoplar K272R/N453D; pTrcPoplar K272R/N453D/C497W; and pTrcPoplar 272R/497W. These plasmids were transformed into BL21(λDE3)pLysS chemically competent cells (Novagen) for analysis. The variants were analyzed for headspace activity (production of isoprene from whole cells), solubility, and specific activity.

II. Cell Growth and Isoprene Production

[0336] The variants, the parent Poplar strain, and the strain containing pTrcKudzu were grown overnight at 37°C in 5 ml of Luria Bertani medium containing either carbenicillin (50 µg/ml-Poplar strains) or kanamycin (50 µg/ml-Kudzu strain). These cultures were diluted into TM3 broth to an OD₆₀₀ of 0.05, supplemented with 10g/L glucose and either 50 µg/ml carbenicillin (Poplar mutants and wild type) or 50 mg/L kanamycin (BL21/pTrcKudzu). The recipe for TM3 broth is as follows: K₂HPO₄ (13.6 g/l) KH₂PO₄ (13.6 g/l), MgSO₄ * 7H₂O (2 g/l) Citric Acid Monohydrate (2 g/L) Ferric Ammonium Citrate (0.3 g/L) (NH₄)₂SO₄ (3.2 g/L) yeast extract (0.2 g/L) 1 ml of 1000x Trace Elements solution, pH adjusted to 6.8 with ammonium hydroxide qs to volume with sterile diH₂O and filter sterilized with a 0.22 micron filter. The recipe for 1000X Trace Elements solution is as follows: Citric Acids * H₂O (40 g/L), MnSO₄ * H₂O (30 g/L), NaCl (10 g/L), FeSO₄ * 7 H₂O (1 g/L), CoCl₂ * 6 H₂O (1 g/L), ZnSO₄ * 7 H₂O (1 g/L), CuSO₄ * 5 H₂O (100 mg/L), H₃BO₃ (100 mg/L), NaMoO₄ * 2 H₂O (100 mg/L). Each component was dissolved one at a time in diH₂O, pH adjusted to 3.0 with HCl/NaOH, qs to volume and filter sterilized with a 0.22 micron filter.

[0337] The diluted culture volume was 25 ml in a 250 ml baffled Bellco Delong flask for growth at 30°C with agitation (225 rpm). Samples were taken aseptically, as indicated, for optical density measurements at A₆₀₀. Two sets of cultures were set up, one for induction with 0.2 mM IPTG and one that remained un-induced. After 3 hours of growth at 30°C with shaking at 200 rpm (OD₆₀₀~0.5), one set of the cultures was induced with 0.2 mM IPTG and incubated for a further 3 h at 30°C with shaking at 200 rpm, the un-induced set was further incubated for the same amount of time. The OD₆₀₀ was determined for all cultures prior to the induction time (3 h post inoculation) and at the time of the measurement of isoprene by the Headspace assay (3 h post-induction, total of 6 h of growth). The cell cultures were centrifuged at 7000 x g for 15 minutes in a Sorvall superspeed centrifuge to pellet the cells. The supernatant was removed and the cell pellet frozen for use in an in vitro assay for isoprene synthase activity. Results of the growth and headspace assays are shown in the following tables.

Table 6-1. Growth of *E. coli* Strains Expressing Poplar IspS Variants
(OD₆₀₀ values)

Strain	Pre-induction (after 3h growth)	Post-induction (3h with IPTG)	Non-induced (after 6h growth)
pTrcPoplar	0.49	3.3	4.7

K272R			
pTrcPoplar K272R/N453D	0.48	3.7	4.8
pTrcPoplar K272R/N453D/C497W	0.41	3.1	5.1
pTrcPoplar K272R/497W H2	0.45	0.8	0.97
pTrcPoplar WT	0.49	3.5	4.7
pTrcKudzu WT	0.44	3.3	4.2

Table 6-2. Production of Isoprene by *E. coli* Strains Expressing Poplar IspS Variants ($\mu\text{g/L}$)

Strain	Pre-induction (after 3h growth)	Post-induction (3h with IPTG)	Non-induced (after 6h growth)
pTrcPoplar K272R	n/a	0	0
pTrcPoplar K272R/N453D	n/a	0	0
pTrcPoplar K272R/N453D/C497W	n/a	0	0
pTrcPoplar K272R/497W H2	n/a	0	0
pTrcPoplar WT	n/a	0	0
pTrcKudzu WT	n/a	1.3	0.29

[0338] Surprisingly, the strain expressing the kudzu IspS demonstrated any measurable isoprene production. This is unexpected given that kinetic properties of the poplar enzymes were reported in the literature to be superior to those of kudzu enzymes. In particular, the prior art describes the specific activity (U/mg) and K_m (μM) of recombinant kudzu IspS to be 0.075 and 7,700 respectively, native aspen IspS to be 0.5 and 8,000 respectively, and recombinant poplar IspS to be 0.16 and 9,000 respectively (Silver and Fall, J Biol Chem, 270:13010-1316, 1995; Miller et al., Planta, 213:483-487, 2001; and Sharkey et al., Plant Physiology, 137:700-712, 2005). The published K_m values for the three enzymes are all quite high and within range of each other, but the specific activity for kudzu isoprene synthase is significantly worse than that of the other two isoprene synthases.

III. Assays for Isoprene Synthase Activity and Solubility

[0339] By using the DMAPP assay, the activity of isoprene synthase can be measured directly as DMAPP is the direct substrate for the enzyme. The cell pellets of the poplar parent and mutant strains, as well as the wild type kudzu were thawed and resuspended in 2 ml PEB (50 mM Tris-HCl, pH 8.0, 20 mM MgCl₂, 2 mM dithiothreitol, and 50% [v/v] glycerol). Cells were lysed by French pressure cell disruption, one pass, at 20,000 psi. The lysate (1 ml) was then centrifuged in a microfuge for 20 min at 20,000 rpm at 4°C. The supernatant was removed and the pellet resuspended in 1 ml of PEB. The supernatant and pellet samples were analyzed by SDS-PAGE, and DMAPP assay, while the total protein content was determined by BCA.

Table 6-3. DMAPP Assay of Isoprene Production from the Supernatant of the Centrifuged Cell Lysate

Strain	OD600 (prior to lysis)	Total Protein (mg/ml)	Isoprene/ Total Protein
Induced			
pTrcPoplar K272R	3.3	1.70	0.03
pTrcPoplar K272R/N453D	3.6	1.3	0.08
pTrcPoplar K272R/N453D/C497W	3.1	1.10	0.08
pTrcPoplar K272R/497W H2	0.803	0.80	0.11
pTrcPoplar WT	3.5	1.30	0.06
pTrcKudzu WT	3.3	1.50	11.15
Uninduced			
pTrcPoplar K272R	4.7	1.7	0.26
pTrcPoplar K272R/N453D	4.8	1.8	0.07
pTrcPoplar K272R/N453D/C497W	5.2	1.9	0.02
pTrcPoplar K272R/497W H2	0.969	1.0	0.17
pTrcPoplar WT	4.6	1.7	0.05
pTrcKudzu WT	4.2	1.9	1.62

[0340] Production was normalized to total cell lysate supernatant protein.

Table 6-4. DMAPP Assay of Isoprene Production from the Pellet of the Centrifuged Cell Lysate

Strain	OD600 (prior to lysis)	Total Protein (mg/ml)	Isoprene/ Total Protein
Induced			
pTrcPoplar K272R	3.3	1.02	0.118
pTrcPoplar K272R/N453D	3.6	1.36	0.000
pTrcPoplar K272R/N453D/C497W	3.1	1.49	0.040
pTrcPoplar K272R/C497W H2	0.80	1.38	0.043
pTrcPoplar WT	3.5	1.57	0.050
pTrcKudzu WT	3.3	1.47	0.040
Uninduced			
pTrcPoplar K272R	4.7	1.40	0.170
pTrcPoplar K272R/N453D	4.8	1.53	0.120
pTrcPoplar K272R/N453D/C497W	5.2	1.42	0.131
pTrcPoplar K272R/C497W H2	0.97	1.55	0.080
pTrcPoplar WT	4.6	1.38	0.120
pTrcKudzu WT	4.2	1.55	0.120

[0341] Production was normalized to total cell lysate pellet protein.

[0342] The poplar variant K272R/C497W showed a 1.8x increase in activity as compared to the wild type in the supernatant fraction of the induced cultures. Likewise, the poplar variants K272R and K272R/C497W showed a 5.2x and 3.4x increase in activity as compared to the wild type in the supernatant fraction of the uninduced cultures. Moreover the poplar variant K272R showed a 2x increase in activity as compared to wild type in the pellet of the induced cultures. However, the most striking result was that the kudzu IspS is more active than the poplar IspS employed herein (185x). In the above tables, H2 is the name of the clone designated pTrcPoplar K272R/C497W H2.

EXAMPLE 7

Subcloning of Kudzu Isoprene Synthase

[0343] In this Example, methods used in the construction of kudzu isoprene synthase (IspS) SELs are described. To create an expression vector for construction of site evaluation libraries (SEL), the kudzu isoprene synthase gene was subcloned into the pET24d vector (Novagen) from the pCR2.1 vector (Invitrogen). The kudzu IspS gene was amplified from pTrcKudzu template DNA using primers MCM50 5'-GATCATGCAT TCGCCCTTAG GAGGTAAAAA AACATGTGTG CGACCTCTTC TCAATTTACT (SEQ ID NO:20); and MCM53 5'-CGGTCGACGG ATCCCTGCAG TTAGACATAC ATCAGCTG (SEQ ID NO:21). PCR reactions were carried out using Taq DNA Polymerase (Invitrogen), and the resulting PCR product was cloned into pCR2.1-TOPO TA cloning vector (Invitrogen), and transformed into *E. coli* Top10 chemically competent cells (Invitrogen). Transformants were plated on L agar containing carbenicillin (50 µg/ml) and incubated overnight at 37°C. Five ml Luria Broth cultures containing carbenicillin 50 µg/ml were inoculated with single transformants and grown overnight at 37°C. Five colonies were screened for the correct insert by sequencing of plasmid DNA isolated from 1 ml of liquid culture (Luria Broth) and purified using the QIAprep Spin Mini-prep Kit (Qiagen). The resulting plasmid, designated MCM93, contains the kudzu IspS coding sequence in a pCR2.1 backbone (Figure 9). The sequence of MCM93 (SEQ ID NO:22) is shown in Figure 10.

[0344] The kudzu coding sequence was removed by restriction endonuclease digestion with *PciI* and *BamHI* (Roche) and gel purified using the QIAquick Gel Extraction kit (Qiagen). The pET24d vector DNA was digested with *NcoI* and *BamHI* (Roche), treated with shrimp alkaline phosphatase (Roche), and purified using the QIAprep Spin Mini-prep Kit (Qiagen). The kudzu IspS fragment was ligated to the *NcoI/BamHI* digested pET24d using the Rapid DNA Ligation Kit (Roche) at a 5:1 fragment to vector ratio in a total volume of 20 µl. A portion of the ligation mixture (5 µl) was transformed into *E. coli* Top 10 chemically competent cells and plated on L agar containing kanamycin (50 µg/ml). The correct transformant was confirmed by sequencing and transformed into chemically competent BL21(λDE3)pLysS cells (Novagen). A single colony was selected after overnight growth at 37°C on L agar containing kanamycin (50 µg/ml). A map of the resulting plasmid designated as pET24D-Kudzu is shown in Figure 11. The sequence of pET24D-Kudzu (SEQ ID NO:23) is shown in Figure 12. IspS Activity was confirmed using the headspace assay as described in Example 5).

EXAMPLE 8

Construction of Isoprene Synthase Site Evaluation Libraries (SELs)

[0345] In order to improve the kinetic parameters of a plant IspS SELs are prepared at sites selected from homology models of both the kudzu and the poplar IspS enzymes. While it is predicted from the homology models that engineering at the indicated sites would improve both enzymes, in this embodiment, kudzu SELs are described. Kudzu IspS surface sites of interest include but are not limited to: 26 L, 30 E, 31 F, 33 Q, 35 L, 36 E, 37 N, 39 L, 40 K, 41 V, 43 K, 44 L, 61 R, 62 V, 63 D, 65 Q, 87 K, 94 E, 95 N, 99 L, 100 D, 105 N, 137 K, 138 E, 143 G, 144 E, 182 N, 184 L, 185 K, 187 G, 189 N, 190 T, 225 P, 226 H, 247 K, 257 T, 258 E, 259 M, 266 D, 334 N, 353 D, 357 S, 358 I, 361 E, 389 N, 392 I, 393 I, 398 K, 401 E, 421 C, 423 Q, 424 Q, 425 E, 426 D, 430 H, 432 L, 433 R, 434 S, 437 D, 443 R, 462 L, 463 E, 476 H, 478 N, 479 D, 485 Q, 508 D, 513 P, 515 A, 532 Q, 533 Y, 537 L, 538 G, 539 R, 542 Y, 543 A, and 557 P. Kudzu IspS active site positions of interest include but are not limited to: 24 P, 25 N, 309 Y, 310 D, 377 L, 381 F, 384 E, 399 Y, 402 N, 403 A, 406 S, 407 S, 409 G, 411 A, 413 L, 449 F, 456 A, 457 T, 458 S, 458 S, 459 A, 460 A, 461 E, 462 L, 463 E, 464 R, 465 G, 466 E, 467 T, 468 T, 469 N, 523 M, 527 S, and 531 Y. Additional kudzu IspS active site positions of interest include but are not limited to: 20 A, 21 N, 22 Y, 23 Q, 271 R, 278 W, 299 F, 302 V, and 408 S. Each library SEL contains clones, maximally including 20 different variants. For example, kudzu isoprene synthase SEL 531 contains variants in which the DNA triplet coding for tyrosine at position 531 of the mature kudzu enzyme is replaced by another DNA triplet encoding: alanine, aspartic acid, cysteine, glutamic acid, phenylalanine, glycine, histidine, isoleucine, lysine, leucine, methionine, asparagine, proline, glutamine, arginine, serine, threonine, valine, or tryptophan. Briefly, DNA triplets of specific positions in the DNA coding strand of the mature IspS are replaced. The mutated IspS nucleic acids are subsequently ligated to a suitable expression vector and used to transform suitable host cells.

[0346] Site evaluation libraries are created either by ordering synthetic constructs (e.g., DNA2.0) or by ordering primers with the “nns” sequence in place of the codon to be mutated. The primers are then be used to mutate the gene to produce an SEL at the indicated site using commercially available mutagenesis kits (e.g., Stratagene) as has been described (e.g., WO0507682A2). The mutated codons are identified by sequence analysis. The site libraries are arrayed in 96 well master plates, and frozen for later use. Cultures are grown from the master plates and prepared for screening.

[0347] The desired end products are IspS enzymes that function optimally in a host metabolically engineered to maximize carbon flow through IspS. To this end several stages of

screening are used to ensure that correct parameters are being addressed. Exemplary screens include but are not limited to: expression, DMAPP feeding for production of HG, microreactor, protein determination, and headspace assays. Expression screen: One example of a method to analyze the level of protein expression is as follows. Soluble and insoluble fractions of cell lysates (obtained from lysed cell cultures) are prepared by centrifugation. The resulting supernatants and pellets are analyzed by SDS-PAGE. The percent soluble protein is determined by densitometry analysis of the protein present in the supernatant versus the pellet.

[0348] In an exemplary embodiment, kudzu site evaluation libraries are constructed in the pET24D expression vector. The pET24D-Kudzu vector, containing the kudzu isoprene synthase gene, serves as the template DNA.

Materials:

pET24D-Kudzu vector (~50 ng/μl)

Kudzu IS site-directed mutagenic primers (Integrated DNA Technologies)

QUIKCHANGE® Multi Site-Directed Mutagenesis Kit (Stratagene)

MJ Research PTC-200 Peltier Thermal Cycler (Bio-Rad Laboratories)

One Shot TOP10 competent cells (Invitrogen)

QIAprep Spin Miniprep Kit (Qiagen)

BL21(λDE3) pLysS competent cells (Invitrogen)

Luria Broth (LB) agar plates

Methods:

[0349] The method of mutagenesis was based on the codon-specific mutation approach, in which the creation of all possible mutations in a specific DNA triplet was performed using a single forward primer with a length of 25 to 45 nucleotides, enclosing a specific designed triple DNA sequence NNS (N = A, C, T or G and S = C or G) corresponding with the sequence of the codon to be mutated. This method results in the random incorporation of nucleotides at a specific pET24D-kudzu codon of interest. Table 8-1 lists the oligonucleotide primers used for mutagenesis, with the number in the primer name corresponding with the codon position in the mature kudzu isoprene synthase enzyme sequence. All oligonucleotide primers were synthesized (Integrated DNA Technologies) on a 100 nmole scale and PAGE purified.

Table 8-1. Kudzu IspS Codon-Specific Mutation Primers

Name	SEQ ID	Primer Sequence
IS_A20	NO:24	CATAATTCCTCGTTCNNNSAACTATCAGCCAAACCTG
IS_N21	NO:25	CATAATTCCTCGTTCGCGANNSTATCAGCCAAACCTGTG
IS_Y22	NO:26	CCCGTCGTTCGCAAACNNSCAGCCAAACCTGTGGAATTC
IS_Q23	NO:27	GTCGTTCGCAAACCTATNNSCCAAACCTGTGGAATTC
IS_R271	NO:28	CTGGATTTTGTACGCGACNNSTGATGGAAGTTTATTC
IS_W278	NO:29	CTGATGGAAGTTTATTCNNSGCACTGGGTATGGCGCC
IS_F299	NO:30	CAAAGCTGTTACTAAAATGNNSGGTCTGGTGACGATCATC
IS_V302	NO:31	CTAAAATGTTTGGTCTGNNSACGATCATCGATGACGTG
IS_S408	NO:32	GAAAACGCCAGCGTTTCCTCCNNSGGTGTAGCGCTGCTGGC

[0350] A PCR reaction was set up in a 0.5ml thin-walled PCR tube following the manufacturer's protocol for the QUIKCHANGE® Multi Site-Directed Mutagenesis Kit (Stratagene): 1µl pET24 Kudzu vector (50ng/µl); 1µl Kudzu IS site-directed forward mutagenic primer (10µM); 2.5µl 10X QUIKCHANGE® Multi Reaction buffer; 1µl dNTP Mix, 1µl QUIKCHANGE® Multi enzyme blend (2.5U/µl); and 18.5µl distilled autoclaved water to provide a 25µl total reaction mix. The pET24 Kudzu SELs were amplified using the following conditions: 95°C, for 1 min (1st cycle only), followed by 95°C for 1 min, 55°C for 1 min, 65°C for 12 min, and repeat cycling 29 times. Then the reaction mixture was subjected to *DpnI* digestion (supplied with QUIKCHANGE® Multi Site-Directed Mutagenesis Kit) by addition of 1.5µl *DpnI* restriction enzyme to each tube, and incubated at 37°C for 2 hours to digest the parental pET24D-kudzu vector. The *DpnI*-treated PCR reaction was then transformed into One Shot TOP10 competent cells (Invitrogen), plated onto LB agar plates containing 50µg/ml kanamycin, and incubated overnight at 37°C. The next day, 96 random colonies were picked and sequenced to identify a minimum of 15 of the possible 19 amino acid variants. Upon identification of the site-directed variants, each variant clone was then inoculated in a 5ml tube of LB + 50µg/ml kanamycin and grown overnight at 37°C with shaking (250 rpm). The following day plasmid DNA was purified using the QIAprep Spin Miniprep Kit (Qiagen). The variants were then transformed into One Shot BL21(λDE3) pLysS competent cells (Invitrogen) for protein expression screening, plated on LB agar plates containing 50µg/ml kanamycin and 30µg/ml chloramphenicol and incubated overnight at 37°C.

[0351] An alternative method for producing pET24D-Kudzu SELs in *E. coli* BL21(λ DE3) pLysS cells was also successfully employed. The TOP10 competent cell transformants obtained from the *DpnI*-treated PCR reaction described above were harvested by applying 3ml of LB media to the top of the agar and resuspending the cells by scraping with a sterile plate spreader. The 3ml of pooled, resuspended cells were then used to inoculate a 25ml shake flask containing LB + 50 μ g/ml kanamycin. The pooled culture was then grown overnight at 37°C with shaking (250rpm). The following day plasmid DNA was purified from the pooled cultures using the QIAprep Spin Miniprep Kit (Qiagen). The pooled plasmid DNA was then transformed into One Shot BL21(λ DE3) pLysS competent cells for protein expression screening as described above.

[0352] To make a master plate, the correct constructs are arrayed in quadruplicate in 96 well plates. One colony of the correct sequence is used to inoculate 4 wells and the plates are grown for several hours to overnight at 37°C in LB containing 50 μ g/ml kanamycin with shaking (200 rpm). Sterile glycerol is added to the cultures to a final concentration of 15% (for a final total volume of 150-200 μ l/well). The plates are then sealed using BREATHE-EASIER (EMS Catalog No. 70536-20) membranes and stored at -80°C.

EXAMPLE 9

Production and Purification of Isoprene Synthase Inclusion Bodies

[0353] Inclusion bodies containing kudzu isoprene synthase were formed when the enzyme is overexpressed in the presence of the chaperone GroELS in the strain BL21(λ DE3). Briefly pETNHISKudzu (U.S. Application No. 61/013,574, herein incorporated by reference) was subcloned into pGro7 (Takara Catalog No. 3340) according to the manufacturer's instructions. A 500 mL of culture was grown essentially as described (Whittington et al., Proc Natl Acad Sci USA, 99:15375-15380, 2002). Despite the presence of chaperone and low temperature of cultivation the culture yielded predominantly inclusion bodies and only low levels of soluble active protein. The inclusion bodies were harvested using the IFOLD Protein Refolding System (Novagen Catalog No. 71552-3) according to the manufacturer's instructions. This procedure led to a high yield (>50 mg) of recombinant kudzu isoprene synthase. The purity of the inclusion body is shown in Figure 13. This preparation was used for the production of rabbit polyclonal anti-isoprene synthase antisera (Invitrogen).

EXAMPLE 10

High Throughput Biochemical Screen of Isoprene Synthase Variants

[0354] This example describes high throughput methods for the determination of isoprene synthase activity. Libraries of BL21(λ DE3)pLysS *E.coli* host cells capable of expressing isoprene synthase variants are arrayed in 96-well plates and stored frozen at -80°C as 15% glycerol stocks as described above in Example 8. To analyze a plate of up to 96 variants, a replica stamp of the glycerol stock master plate is made with a 96-pin MULTI-BLOT floating pin tool (V&P Scientific Catalog No. VP 408AF) onto Luria broth agar containing appropriate antibiotic(s) (e.g., 30 $\mu\text{g/mL}$ chloramphenicol, 50 $\mu\text{g/mL}$ kanamycin). The replica plate is incubated over night at 30°C to allow growth of bacterial patches. Using the same floating pin replicator a 96-square deep well plate containing 250 μL of TM3 medium supplemented with 0.08% Biospringer yeast extract and 1% glucose plus antibiotics (30 $\mu\text{g/mL}$ chloramphenicol, 50 $\mu\text{g/mL}$ kanamycin) is inoculated from the agar plate and incubated overnight at 30°C . The recipe for TM3 broth is as follows: K_2HPO_4 (13.6 g/l) KH_2PO_4 (13.6 g/l), $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ (2 g/l) Citric Acid Monohydrate (2 g/L) Ferric Ammonium Citrate (0.3 g/L) $(\text{NH}_4)_2\text{SO}_4$ (3.2 g/L) yeast extract (0.2 g/L) 1 ml of 1000x Trace Elements solution, pH adjusted to 6.8 with ammonium hydroxide qs to volume with sterile dH_2O and filter sterilized with a 0.22 micron filter. The recipe for 1000X Trace Elements solution is as follows: Citric Acids $\cdot \text{H}_2\text{O}$ (40 g/L), $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ (30 g/L), NaCl (10 g/L), $\text{FeSO}_4 \cdot 7 \text{H}_2\text{O}$ (1 g/L), $\text{CoCl}_2 \cdot 6 \text{H}_2\text{O}$ (1 g/L), $\text{ZnSO}_4 \cdot 7 \text{H}_2\text{O}$ (1 g/L), $\text{CuSO}_4 \cdot 5 \text{H}_2\text{O}$ (100 mg/L), H_3BO_3 (100 mg/L), $\text{NaMoO}_4 \cdot 2 \text{H}_2\text{O}$ (100 mg/L). Each component was dissolved one at a time in dH_2O , pH adjusted to 3.0 with HCl/NaOH, qs to volume and filter sterilized with a 0.22 micron filter. The overnight cultures are diluted with the same medium to an OD_{600} of 0.05 and grown in another 96-square deep well plate (Thomson Instrument, Catalog No. 951652C), with each well containing 600 μL of the dilution. The dilutions are grown at 30°C with shaking to an OD_{600} of 0.8 and are then induced with IPTG added to a concentration of 400 μM . The plate is grown for 5 hours and OD_{600} is determined for quality control and normalization.

[0355] A volume of 400 μL of culture is transferred into a new 96-well plate (Perkin Elmer, Catalog No. 6008290) and cells are harvested by centrifugation in a Beckman Coulter Allegra 6R centrifuge at 2500 x g. The pellet is resuspended in 200 μL of hypotonic buffer (5mM MgCl_2 , 5mM Tris HCl, 5mM DTT pH 8.0) and the plate is frozen at -80°C for a minimum time of 60 min. Cell lysate is prepared by thawing the plate and adding 32 μL of isoprene synthase DMAPP assay buffer (57mM Tris HCl, 19 mM MgCl_2 , 74 $\mu\text{g/mL}$ DNase I

(Sigma Catalog No. DN-25), 2.63×10^5 U/mL of READYLYSE lysozyme solution (Epicentre Catalog No. R1802M), and 5 mg/mL of molecular biology grade BSA. The plate is incubated with shaking at 25°C for 30 min and then placed on ice. For isoprene production, an 80 μ L aliquot of lysate is transferred to a 96-deep well glass plate (Zinsser Catalog No. 3600600) and 20 μ L of a 10 mM DMAPP solution in 100mM KHPO₄, pH 8.2 (Cayman Chemical Catalog No. 63180) is added. The plate is sealed with an aluminum plate seal (Beckman Coulter Catalog No. 538619) and incubated with shaking at 30°C of 60 minutes. The enzymatic reactions are terminated by heating the glass block (70°C for 5 min). The headspace of each well is quantitatively analyzed as described in Example 5.

[0356] To determine protein concentration 5 μ L or more of lysate is run on precast gels (Invitrogen Catalog No. NP0301BOX) for western blot analysis after transfer to a nitrocellulose membrane (Invitrogen Catalog No. LC2000). The primary antibody employed is an anti-isoprene synthase antibody of Example 9. Primary antibody binding is followed by development with a secondary antibody labeled with Alexa Fluor 488 (Invitrogen Catalog No. A-11008) to permit quantitative signal determination. The western blot procedure was carried out as described by Invitrogen. The fluorescence signal was recorded with a Molecular Dynamics STORM instrument using the blue filter setting and quantitatively analyzed with the Molecular Dynamics IMAGEQUANT image analysis software package. Specific activity of the library members was calculated from the ratio of the amount of isoprene produced divided by either the A₆₀₀ of the induction cultures or the isoprene synthase protein concentration determined by western blot. Isoprene synthase protein standard was calibrated by standard gel densitometry with BSA stained with Coomassie brilliant blue R250 serving as primary standard. Increased, decreased, or no change in specific activity of the entire library was tabulated for further analysis. Figure 14 provides graphs showing isoprene synthase activity of kudzu site evaluation library (SEL) members for positions Y22, A20, and S408. Most members show highly decreased activity relative to wild type, while conservative substitutions show a lesser decrease in activity. . Activity of variant A20G approximates that of the wild type kudzu enzyme, indicating that it is a candidate partner for a combinatorial mutant. Interestingly, variant S408D of library S408 showed an increase in activity compared to wild type thus providing another candidate partner for a combinatorial mutant.

Example 11

Isoprene Synthase Truncations

[0357] This example describes the identification of the amino acid sequence of the protein in the lower band of the doublet seen in purified poplar IspS preparations (see Figure 21). A series of N-terminally truncated IspS molecules based on putative cleavage sites identified by mass spectrometry was also generated. A shorter N-terminal truncation of IspS (the “MEA” truncation in pDu39, see below) was also generated, to examine the effect of further truncation on IspS activity (Williams DC, McGarvey DJ, Katahira EJ, Croteau R (1998) *Biochemistry* 37:12213–12220).

I. Construction of an N-terminally 6XHis-tagged IspS (in pDu27) for protein purification:

[0358] The full length *P. alba* IspS from the template *P. alba* pET24a (Figures 19 and 20) was prepared by PCR. The following PCR reaction was prepared: 1 µl (Template)-*P. alba* pET24a, 5 µl 10X PfuUltraII Fusion buffer, 1 µl dNTP's (10 mM), 1 µl primer (50 µM) primer F-(MCM219), 1 µl primer (50 µM) primer R-(MCM182), 41 µl diH₂O and 1 µl of PfuUltra II Fusion DNA Polymerase (Stratagene). PCR cycling parameters were as follows: 95 °C 1 min., 95 °C 1 min, 55 °C 20 sec., 72 °C 27 sec. for 29 cycles followed by 72 °C 3 min and 4 °C until cool, using an Eppendorf Mastercycler. The PCR product was gel extracted and purified, using 0.8% E-gel (Invitrogen) and Qiagen QIAquick Gel Extraction and QIAprep Spin Miniprep kits, according to the manufacturer's recommended protocol. A 3 µl aliquot of purified product was ligated to the pET200D/TOPO vector (Invitrogen), according to the manufacturer's protocol. The reaction was incubated for 5 minutes at room temperature, and the 6 µl topoisomerase mixture was then transformed into *E. coli* Top10 chemically competent cells (Invitrogen) according to the manufacturer's protocol. Transformants were selected for on LB plates containing kanamycin (50 µg/ml) (Kan50), and incubated at 37 °C overnight. Five colonies were picked and screened using PuReTaq Ready-To-Go PCR Beads (Amersham) using the T7 Forward and MCM182 primers. Clones harboring inserts of the correct size were further verified by sequencing using the T7 Forward and T7 Reverse primers (Quintara Biosciences). One construct, pDu27 (Figures 16-18), was chosen for further study. A 1 µl aliquot of the plasmid preparation was transformed into BL21(λDE3)*pLysS* (Invitrogen) according to the manufacturer's protocol. Transformants were selected for on LB plates containing Kan50+ and chloramphenicol (35 µg/ml) (Cm35) and incubated at 37 °C overnight. The resulting strain was used for expression and purification of N-terminally 6XHis-tagged *P. alba* IspS.

[0359]

II. Purification of 6XHis-tagged IspS

Expression of 6XHis-tagged IspS

[0360] N-terminally 6XHis-tagged IspS was expressed and purified from strain MD08-99. The growth procedure is suitable for histidine tagged enzymes expressed in BL21(λ DE3)*pLysS* cells. A 10 ml of overnight culture was prepared for each 1L of planned growth. The appropriate antibiotics (50 mg/ml kanamycin, 50 mg/ml chloramphenicol, and/or 50 mg/ml Carbenecillin) was added to 10 ml of LB medium in a 25 ml flask and was inoculated with 1 colony from a fresh plate of cells or directly from glycerol frozen cell stock. Cultures were grown at 30 °C overnight with shaking at ~220 rpm. Day cultures were prepared in 1 liter of LB medium with appropriate antibiotics for each culture. Each 1L day culture was inoculated with 10 ml of overnight culture and grown at 30-37 °C with shaking at ~220 rpm until the OD₆₀₀ reached ~0.4-0.6. Day cultures were then induced with 400 μ M IPTG and allowed to continue growing at 30°C with shaking at 220 rpm for ~5-6 hours. Cells were then harvested by centrifugation at 10,000 x g for 10 min, 4 °C. Following Harvest, cells were used directly or stored at -80°C until ready to process.

Purification of 6XHis-tagged IspS

[0361] For purification of histidine tagged enzymes from BL21(λ DE3)*pLysS* cells, cells were gently resuspended in fresh Lysis buffer (Lysis buffer: Ni wash buffer + 0.5 mM PMSF, 0.01% Tween-20, 1 mg/ml lysozyme, 0.2 mg/ml DNaseI; Ni wash buffer: 50 mM NaH₂PO₄, 300 mM NaCl, 20 mM Imidazole, pH 8.0). Approximately 40-50 ml of lysis buffer was used per 1L of cell pellet. Cells were then incubated on ice for approximately 30 min. The cell suspension was then lysed fully by passing 2-3 times through a french pressure cell (large french press cell at 1200 psi/High setting) until lysate started to look clear. A sample of the lysate was saved for activity assay and gel analysis (~100 μ l). The lysate was then clarified by centrifuging the lysate at 30,000xg for 30 min, 4°C in a Sorvall Discovery 90SE ultracentrifuge. The supernatant was removed and retained. A sample of the “clarified lysate” was saved for activity assay and gel analysis (~100 μ l).

[0362] The clarified lysate was run over HisTrap HP columns (GE Healthcare) using a gradient from 0-100% Ni buffer B. Samples were then analyzed by SDS-PAGE gel (4-12% gel NUPAGE, Invitrogen) according to manufacturer’s directions. Desired fractions were concentrated on spin filters (Vivaspin-20, Sartoris,) and then desalted over a HiPrep 26/10

Desalting column (GE heathcare) packed with Sephadex G25 resin. The G-25 buffer consisted of 50 mM HEPES, 50 mM NaCl, and 1 mM DTT, pH 7.4. The desired sample was then purified over a HiTrap Q HP column (GE) using a gradient elution from 0% Q seph buffer A to 100% Q seph buffer B (Q seph buffer A: 50 mM Tris, 0.05 M NaCl, 1 mM DTT, pH 7.6 and Q seph buffer B: 50 mM Tris, 1.0 M NaCl, 1 mM DTT, pH 7.6). Fractions containing the desired protein were analyzed and concentrated. Sample buffer was then exchanged into 50 mM HEPES, 50 mM NaCl, pH 7.4 with 1 mM DTT by passing the sample over a Hi Prep 26/10 Desalting column (GE heathcare) packed with Sephadex G25 resin. A final polishing step of Gel filtration was used when necessary. The sample was passed through a Hi Load 26/60 Superdex 200 prep grade (GE) in gel filtration buffer: (50 mM HEPES, 150 mM NaCl, 1mM DTT, pH 7.4). Fractions were then analyzed and concentrated. The samples were then stored at -80 °C. For preparation for analysis of the band, the sample is run on an SDS-PAGE gel (4-12% NUPAGE gel, Invitrogen), stained and the desired band excised and processed as described below.

III. Mass Spectrometry of Isoprene Synthase

Sample preparation

[0363] An In-Gel Digestion and LCQ-Deca Mass Spectrometry Procedure was utilized (Modified Rosenfeld in-gel Digest Protocol) (Rosenfeld *et al*, Anal. Biochem. (1992) 203, 173-179; Hellman *et al*, Anal Biochem, (1995) 224, 451-455). The purified sample of Isoprene synthase was run on a 4-12% SDS-PAGE (NUPAGE, Invitrogen) and stained with Coomassie Brilliant Blue R-250 (Thermo Scientific) or SimplyBlue Safe Stain (Invitrogen). Band(s) of interest were excised from the gel and destained. Each gel slice was diced into small pieces ~1 mm x 1mm and placed into 0.65 mL “slick” (siliconized) tubes from PGC Scientific. Approximately 100µL of 25mM NH₄HCO₃/50% ACN/H₂O was added to each tube and vortexed for 10 min. Supernatants were extracted and discarded. These steps were repeated twice. Then gel pieces were then run in a Savant SpeedVac to dryness (~ 20 to 30 min).

[0364] Samples were then reduced and alkylated. For reduction, 25 µL (or enough to cover) of 10 mM DTT in 25 mM NH₄HCO₃ (prepared fresh) were added to dried gels. Tubes were then vortexed and spun briefly. Reactions were incubated at 50 °C for 1 hour. For alkylation, supernatants were removed and 25 µL or more of 55 mM iodoacetamide (IAA) in 25 mM NH₄HCO₃ were added to the gel. Reaction tubes were vortexed and spun briefly

again. Reactions were allowed in dark for 1 hour at room temperature. Supernatants were removed and gels were washed with ~ 100 μ L of 25mM NH_4HCO_3 /50% ACN/ H_2O , by vortexing for 10 min and briefly spinning. Supernatant were removed and the wash step was repeated once. Gel pieces were then dried in a SpeedVac (~ 15 -30 min).

[0365] Digestion buffer was prepared by adding 400 μ L of 0.1% n-octyl B-D-glucopyranosidase water to 100 μ L of 8M Urea. 400 μ L of this digestion buffer was added to 20 μ g of freshly prepared Trypsin. 0.05 μ g/ μ L of sequencing-grade Trypsin was prepared from one vial of 20 μ g sequencing grade trypsin (Promega) that was dissolved into 400 μ L of 1.6 M Urea solution. Trypsin enzyme solution was added enough to cover gel pieces. Tubes were covered with parafilm and incubated at 37 °C overnight (16-20 hrs). It was ensured that there is a little extra buffer above the gel.

[0366] Peptides were extracted from gels by briefly vortexing and spinning the digest. The digest solution was transferred with gel loading tips into a 0.65 mL siliconized tube. 50 μ L (enough to cover) of 50% ACN/0.1% FA/ H_2O were added to the gel pieces and samples were vortexed for 10 min, spun, and then sonicated for five min. Extracted peptides were pooled together in one tube. Extraction steps were repeated two to three more times until the gel pieces became white in appearance and shrank in size. Extracted digests were vortexed, spun and dried in a Speed Vac to a volume of 55 μ L. In cases where the volume was less than 55 μ L, enough 0.1% FA was added to make up a final volume of 55 μ L.

Mass Spectrometry

[0367] The sample was injected onto a Thermofinnigan (San Jose, CA) LCQ-Deca electrospray ionization (ESI) ion-trap mass spectrometer. A Vydac C18 column (5 μ , 300A, 0.2 x 150 mm, Michrom Bioresources, Auburn, CA) was used with a flow rate of 200 μ L/min. The injection volume was 50 μ L, and was filtered through an on-line trapping cartridge (Peptide CapTrap, Michrom Bioresources, Auburn, CA) before loading onto the column. Separation of the in-gel digest was performed with the following gradient (Solvent A: 0.1% trifluoroacetic acid in H_2O (J.T. Baker, Phillipsburg, NJ), Solvent B: 0.08% trifluoroacetic acid in acetonitrile (J.T. Baker, Phillipsburg, NJ)):

Table 11-1. Gradient Table

	min	A%	B%

0	0.00	100	0
1	10.00	86	14
2	16.00	81	19
3	20.00	78	22
4	21.00	77	23
5	22.00	75	25
6	24.00	73	27
7	32.00	69	31
8	34.00	66	34
9	37.00	64	36
10	47.00	60	40
11	50.00	30	70
12	55.00	100	0
13	60.00	100	0
14	65.00	100	0
15		100	0

Mass Spectrometry Results

[0368] An aliquot of 6.4 µg of protein was loaded into 5 lanes on a 4-12% bis-tris NUPAGE gel (Invitrogen), MOPS buffer, 50 min run. As described above, the gel was stained for 2 minutes, and then de-stained for 15 minutes. The gel was washed in H₂O, then all bands were excised, cut into small pieces, and destained. Gel pieces were reduced and alkylated with DTT/IAA for 1 hour each at 52 °C and RT, respectively. Trypsin was added for an overnight digestion. Extracted peptides were run on the LCQ-Deca. Figures 21 and 22 show the results for the mass spectrometry analysis. The lower doublet band (in Figure 21) is identified as IspS. Figure 22 shows that N-terminal truncations were observed after amino acids 39, 40, 42, and 44 (according to the peptide sequence of the N-terminally His-tagged IspS protein in pDu27). The C-terminus of IspS is intact.

IV. Construction of N-terminally truncated IspS variants:

[0369] All truncated constructs without affinity tags were generated using the QuickChange Site-directed Mutagenesis kit (Stratagene) using the template P. alba pET24a

for PCR amplification. Approximately 50 ng of template DNA was used for amplification (with an Eppendorf Mastercycler Gradient PCR Machine) of the mutagenized PCR product with the Forward (For) and Reverse (Rev) primer pairs that correspond to each relevant truncation (QC MSV For and QC MSV Rev, for example, see Table 11-2). The following PCR reaction mixtures were used: 1 µl *P. alba* pET24a, 5 µl 10X PfuUltra HF buffer, 1 µl dNTP's, 1 µl (50 µM) primer-For (*e.g.* QC MSV For), 1 µl (50 µM) primer-Rev (*e.g.* QC MSV Rev), 2 µl DMSO, 39 µl diH₂O, 1 µl PfuUltra HF Polymerase (Stratagene). The following PCR cycling parameters were used: 95 °C 1 min, 95 °C 30 sec., 55 °C 1 min., 68 °C 7.3 min. for one cycle followed by 95 °C 30 sec., 55 °C 1 min., 68 °C 7.3 min for a total of 18 cycles and then followed by 4 °C. The PCR products were treated with 1-2 µl of *DpnI* (Roche) for 1-3 hour at 37 °C. A 5 µl aliquot of the *DpnI* treated products was visualized on a 0.8% E-gel (Invitrogen). A 1 µl aliquot of each product was transformed into chemically competent *E. coli* Top10 cells (Invitrogen) according to the manufacturer's protocol. Transformants were selected for on LB medium containing kanamycin at a concentration of 50 µg/ml (Kan50), and incubated overnight at 37 °C. Five colonies of each transformation were selected and grown to stationary phase in 3 ml liquid LB Kan50. Plasmids were purified using a Qiagen miniprep kit according to the manufacturer's recommended protocol. Purified plasmids were sequenced (by Quintara Biosciences) with T7 Forward and Reverse primers, and confirmed for their respective truncation. The resulting plasmids (pDU39 through pDU43, see Table 11-4, Figures 23-34) were transformed into chemically competent *E. coli* BL21(λDE3) *pLysS* (Invitrogen) according to the manufacturer's recommended protocol. Table X describes the strains used for expression of truncated IspS enzymes.

[0370] Constructs with affinity (6XHis) and proteolysis (TEV, Tobacco Etch Virus) tags were generated using *P. alba* pET24a as a template for PCR reactions. PCR reaction mixtures were prepared as follows: 1ul (*P. alba* pET24a), 5ul 10X PfuUltraII Fusion buffer, 1ul dNTP's (10 mM), 1ul primer (50 µM) Alba FL-NdeI-For or Alba TRC (MEA)-NdeI-F, 1ul primer (50uM) Alba FLTRC (+) TEV-R, 41 ul diH₂O and 1ul of PfuUltra II Fusion DNA Polymerase from Stratagene. PCR cycling parameters were as follows: 95 °C 1 min., 95 °C 30 sec., 55 °C 20 sec., 72 °C 25 sec. for one cycle and then repeating 95°C 30 sec., 55 °C 20 sec., 72 °C 25 sec. for an additional 28 cycles, followed by 72°C 3min and then 4 °C. After amplification and verification of the correct molecular weight of the product by visualization on 0.8% E-gel (Invitrogen), PCR products were digested with restriction enzymes *NdeI* and *XhoI* (Roche) for 2 hours at 37 °C, and then gel purified using the Qiaquick Gel Purification system (Qiagen) according to the manufacturer's recommended protocol. 3 µl of purified

product was ligated to pET-24a (Invitrogen) that was digested with *NdeI* and *XhoI* (Roche), gel purified and dephosphorylated (using SAP, shrimp alkaline phosphatase) (Promega) according to the manufacturer's recommended protocols. T4 ligase (NEB) was used for the ligation reaction, which was incubated overnight at 16 °C. The ligation reaction was dialyzed into water for 30 min., and 2 µl of the reaction were used to electroporate MCM331 (see below) competent cells. Cells were allowed to recover at 30 °C for 2 hours, and then selected on Kan50 with 5 mM (R)-(-)-Mevalonolactone (MVA) (Sigma) spread onto the plate. Positive transformants were inoculated into 3 ml of liquid LB Kan50, and plasmids were isolated using the QIAprep Spin miniprep kit (Qiagen). Inserts were verified by restriction digestion using *NdeI* and *XhoI* (Roche) and positive clones were sequenced (Quintara Biosciences) with T7 promoter and T7 terminator sequencing primers. 1 µl of each plasmid (see Table 11-4 for plasmid description and Figures 35-39) was transformed into chemically competent *E. coli* BL21(λDE3) pLysS (Invitrogen) according to the manufacturer's recommended protocol. Transformants were selected on LB Kan50 + Cm35 (Chloramphenicol 35 ug/ml) plates and incubated at 37 °C. See Table 11-5 for a description of all expression strains.

[0371] Strain MCM331 was prepared as follows. A synthetic operon containing mevalonate kinase, mevalonate phosphate kinase, mevalonate pyrophosphate decarboxylase, and the IPP isomerase was integrated into the chromosome of *E. coli*. If desired, expression may be altered by integrating different promoters 5' of the operon.

i) Target vector construction

[0372] The attTn7 site was selected for integration. Regions of homology upstream (attTn7 up) (primers MCM78 and MCM79) and downstream (attTn7 down) (primers MCM88 and MCM89) were amplified by PCR from MG1655 cells. A 50 µL reaction with 1 µL 10 µM primers, 3 µL ddH₂O, 45 µL Invitrogen Platinum PCR Supermix High Fidelity, and a scraped colony of MG1655 was denatured for 2:00 at 94 °C, cycled 25 times (2:00 at 94 °C, 0:30 at 50 °C, and 1:00 at 68 °C), extended for 7:00 at 72 °C, and cooled to 4 °C. This resulting DNA was cloned into pCR2.1 (Invitrogen) according to the manufacturer's instructions, resulting in plasmids MCM278 (attTn7 up) and MCM252 (attTn7 down). The 832 bp *ApaI*-*PvuI* fragment digested and gel purified from MCM252 was cloned into *ApaI*-*PvuI* digested and gel purified plasmid pR6K, creating plasmid MCM276. The 825bp *PstI*-*NotI* fragment digested and gel purified from MCM278 was cloned into *PstI*-*NotI* digested and gel purified MCM276, creating plasmid MCM281.

ii) Cloning of lower pathway and promoter

[0373] MVK-PMK-MVD-IDI genes were amplified from pTrcKKDyIkIS with primers MCM104 and MCM105 using Roche Expand Long PCR System according to the manufacturer's instructions. This product was digested with *NotI* and *ApaI* and cloned into MCM281 which had been digested with *NotI* and *ApaI* and gel purified. Primers MCM120 and MCM127 were used to amplify CMR cassette from the GeneBridges FRT-gb2-Cm-FRT template DNA using Stratagene Pfu Ultra II. A PCR program of denaturing at 95 °C for 4:00, 5 cycles of 95 °C for 0:20, 55 °C for 0:20, 72 °C for 2:00, 25 cycles of 95 °C for 0:20, 58 °C for 0:20, 72 °C for 2:00, 72 °C for 10:00, and then cooling to 4 °C was used with four 50 uL PCR reactions containing 1uL ~10ng/μL template, 1 μL each primer, 1.25 μL 10mM dNTPs, 5 μL 10x buffer, 1 μL enzyme, and 39.75 μL ddH₂O. Reactions were pooled, purified on a Qiagen PCR cleanup column, and used to electroporate water-washed Pir1 cells containing plasmid MCM296. Electroporation was carried out in 2 mM cuvettes at 2.5V and 200 ohms. Electroporation reactions were recovered in LB for 3 hr at 30 °C. Transformant MCM330 was selected on LA with CMP5, Kan50.

iii) Integration into *E. coli* chromosome

[0374] Miniprep DNA (Qiaquick Spin kit) from MCM330 was digested with *SnaBI* and used to electroporate BL21(DE3) (Novagen) or MG1655 containing GeneBridges plasmid pRedET Carb. Cells were grown at 30 °C to ~OD1 then induced with 0.4% L-arabinose at 37 °C for 1.5 hours. These cells were washed three times in 4 C ddH₂O before electroporation with 2 μL of DNA. Integrants were selected on L agar with containing chloramphenicol (5 μg/ml) and subsequently confirmed to not grow on L agar + Kanamycin (50 μg/ml). BL21 integrant MCM331 and MG1655 integrant MCM333 were frozen.

Table 11-2. Primers

MCM219	caccatgcggtttagcgtgtcca (SEQ ID NO:33)
MCM182	gggcccggtttaaactttaactagactctgcagtttagcgttcaaacggcagaa (SEQ ID NO:34)
QC MSV For	gaaggagatatacatatgagcgtgtccaccgaaaatg (SEQ ID NO:35)
QC MSV Rev	cattttcgggtggacacgctcatatgtatatctccttc (SEQ ID NO:36)
QC MVS For	gaaggagatatacatatggtgtccaccgaaaatgtgtc (SEQ ID NO:37)
QC MVS Rev	gacacattttcgggtggacaccatattgtatatctccttc (SEQ ID NO:38)
QC MTE For	gaaggagatatacatatgaccgaaaatgtgtctttcac (SEQ ID NO:39)
QC MTE Rev	gtgaaagacacattttcgggtcatatgtatatctccttc (SEQ ID NO:40)
QC MNV For	gaaggagatatacatatgaatgtgtctttcaccgaaac (SEQ ID NO:41)
QC MNV Rev	gtttcgggtgaaagacacattcatatgtatatctccttc (SEQ ID NO:42)
QC MEA For	gaaggagatatacatatggaagctcgtcgtttctgcg (SEQ ID NO:43)
QC MEA Rev	cgcagaacgacgagcttccatattgtatatctccttc (SEQ ID NO:44)
Alba FL-NdeI-For	gaaggagatatacatatgcggtttagcgtg (SEQ ID NO:45)

Alba FLTRC (+) TEV-R	cccgcgcttactcgaggccctgaaaatacagggttttcgcg ttcaaacggcagaatcggtt (SEQ ID NO:46)
Alba TRC (MEA)- NdeI-F	gaaactgaaacccatatggaagctcgctcgttctgc (SEQ ID NO:47)

Table 11-3. Primers for construction of strain MCM331

MCM78	attTn7 up rev for integration construct	gcatgctcgagcgccgcTTTAAATCAAACATCCTGCCAACTC (SEQ ID NO:48)
MCM79	attTn7 down rev for integration construct	gatcgaagggcgatcgTGTCACAGTCTGGCGAAACCG (SEQ ID NO:49)
MCM88	attTn7 up forw for integration construct	ctgaattctgcagatatcTGTTTTTCCACTCTTCGTTCACTT (SEQ ID NO:50)
MCM89	attTn7 down forw for integration construct	tctagagggcccAAGAAAAATGCCCCGCTTACG (SEQ ID NO:51)
MCM104	GI1.2 promoter – MVK	Gatcgcgccgagcgcccttgacgatgccacatcctgagcaa Ataattcaaccactaattgtgagcggataacacaaggagg Aaacagctatgtcattaccgttcttaacttc (SEQ ID NO:52)
MCM105	aspA terminator – yIDI	Gatcgggccccaagaaaaaaggcacgtcatctgacgtgcc Ttttttatttgtagacgcgttggttatagcattcta (SEQ ID NO:53)
MCM120	Forward of attTn7: attTn7 homology, GB marker homology	aaagtagccgaagatgacggtttgtcacatggagtt ggcaggatggttgattaaaagcAATTAACCCCTCACTA AAGGGCGG (SEQ ID NO:54)
MCM127	Rev complement of 1.2 GI: GB marker homology(extra long), promoter, RBS, ATG	AGAGTGTTCACCAAAAAATAAACCTTTCCCGGTGCAGaa Gttaagaacggtaatgacatagctgtttcctccttggtgtt Atccgctcacaattagtggttgaattatttgctcaggatg tggcatcgtaagggtTAATACGACTCACTATAGGGCTCG (SEQ ID NO:55)

Table 11-4. Plasmids for expression of IspS variants

MD09-161	pET24a-P.alba FL C-Term (+) TEV, His tag/MCM331
MD09-163	pET24a-P.alba TRC (MEA) C-Term (+) TEV, His tag/MCM331
pDu27	Alba-FL-pET200/D-TOPO
pDu39	Mtg pET24a-P.alba-MEA /Top10 (Untagged)
pDu40	Mtg pET24a-P.alba-MNV/Top10 (Untagged)
pDu41	Mtg pET24a-P.alba-MSV/Top10 (Untagged)
pDu42	Mtg pET24a-P.alba-MTE/Top10 (Untagged)
pDu43	Mtg pET24a-P.alba-MVS/Top10 (Untagged)

Table 11-5. Strains for expression of IspS variants

MD08-99	Alba-FL-pET200/D-TOPO (pDu27) in BL21 (λ DE3) pLysS
MD09-165	BL21(λ DE3)pLysS, pET24a-P.alba FL C-Term (+) TEV, His tag
MD09-167	BL21(λ DE3) pLysS, pET24a-P.alba TRC (MEA) C-Term (+) TEV, His tag
MD09-173	BL21(λ DE3)pLysS, pET24a-P.alba (MEA) Untagged (pDu39)
MD09-174	BL21(λ DE3)pLysS, pET24a-P.alba (MNV) Untagged (pDu40)
MD09-175	BL21(λ DE3)pLysS, pET24a-P.alba (MSV) Untagged (pDu41)

MD09-176	BL21(λ DE3)pLysS, pET24a-P.alba (MTE) Untagged (pDu42)
MD09-177	BL21(λ DE3)pLysS, pET24a-P.alba (MVS) Untagged (pDu43)

V. Biochemical Analysis of IspS Truncations

[0375] The relative activity of the various N-terminally truncated IspS enzymes was determined by DMAPP assay. The strains described above were analyzed via DMAPP assay in a 96-well plate. All strains were assayed in quadruplicate. The “Full Length” variant refers to the IspS enzyme expressed in BL21(λ DE3) pLysS (Invitrogen) from the P. alba pET24a plasmid.

Table 11-6. DMAPP Assay of N-terminal Truncations

Variant	Activity Data			Avg OD ₆₀₀	OD Normalized
	Average	Std Dev	%CV		
MD09-173	1125.6	93.7	8.3	5.2	217.6
MD09-174	118.6	8.0	6.7	5.2	22.8
MD09-175	1064.6	71.7	6.7	4.9	219.0
MD09-176	1179.1	64.7	5.5	4.9	238.7
MD09-177	831.7	89.6	10.8	4.9	168.2
Full Length	805.8	65.1	8.1	5.0	161.3

[0376] Results: Table X shows that when normalized for OD₆₀₀, strains MD09-173 (with plasmid pDu39), MD09-175 (pDu41), MD09-176 (pDu42), and MD09-177 (pDu43) all displayed higher DMAPP activity than the “Full Length” wild type IspS enzyme (in strain BL21(λ DE3) pLysS with P. alba pET24a).

VI. Detailed Kinetic Analysis of the “MEA” Truncation in IspS

[0377] The relative specific activity was determined and the kinetics of “N-terminally truncated” *P. alba* isoprene synthases were examined compared to “Full length” *P. alba* isoprene synthases. Four strains expressing four different constructs were used in this analysis: BL21(λ DE3) pLysS with P. alba pET24a; MD09-173; MD09-165; and MD09-167 (Strains described above in detail). These strains express “full length” *P. alba* IspS, “truncated” *P. alba* IspS (the MEA truncation), “full length” C-terminally TEV and His-tagged *P. alba* IspS, and “truncated” C-terminally TEV and His-tagged *P. alba* IspS, respectively. In the experiments described below, “truncated” refers specifically to the MEA variant of *P. alba* IspS.

[0378] All strains were inoculated into LB containing 30 mg/L chloramphenicol(Cm) and 50 mg/L kanamycin and grown overnight in 2 mL culture tubes at 37 °C. The overnight cultures were diluted 1:100 in 25 mL of LB broth containing 30 mg/L chloramphenicol(Cm) and 50 mg/L kanamycin the following morning and grown at 37 °C until OD~0.5. Each strain was grown in triplicate. The cultures were then induced with 400 uM IPTG and incubated at 30 °C for 4 hours. 20 mL of each culture were centrifuged at 3000 x g for 20 min. and the supernatant was discarded. The pellets were frozen at -80 °C overnight. Pellets were resuspended in 2 mL of a buffer containing 100 mM Tris, 100 mM NaCl, 0.25 mg/mL lysozyme and 0.25 mg/mL DNAase, pH 8. Cell suspensions were french pressed at 20,000 psi twice and the lysate was then centrifuged at 14000 x g for 20 minutes to yield cell free extract that was used for kinetic studies and protein concentration determination.

[0379] To measure specific activity, 5 µL of cell free extract from each strain was incubated with 5 mM DMAPP, 50 mM MgCl₂ in a buffer containing 100 mM Tris and 100 mM NaCl (pH 8) to a final volume of 100 µL for 15 min. at 30 °C in gas tight 2 mL vials. Reactions were terminated with the addition of 100 µL of 500 uM EDTA, pH 8. Samples were analyzed by GC-MS to determine the concentration of isoprene in the headspace of the vials.

[0380] To determine k_{cat} and K_M , 5 uL of cell free extract from each strain was incubated with DMAPP at concentrations ranging from 0.625 to 40 mM DMAPP in a buffer containing 100 mM Tris, 50 mM MgCl₂ and 100 mM NaCl (pH 8) to a final volume of 100 µL for 15 min. at 30 °C in gas tight 2 mL vials. Reactions were terminated with the addition of 100 µL of 500 mM EDTA, pH 8. Samples were analyzed by GC-MS to determine the concentration of isoprene in the headspace of the vials. Data were analyzed using Kaleidagraph and fit to following equation for uncompetitive substrate inhibition: $rate/E = k_{cat} * S / (K_M + S * (1 + S/K_i))$. All data were run in triplicate with the exception of MD09-167 with 2.5 mM DMAPP which was run in duplicate.

[0381] Cell free extract was run on a Caliper microfluidic electrophoresis instrument (Caliper Life Sciences, Hopkinton, MA, USA) in order to quantify the amount of isoprene synthase in each sample. The microfluidic chip and protein samples were prepared according to the manufacturer's instructions (LabChip® HT Protein Express, P/N 760301). Culture lysates were prepared in 96-well microtiter plates by adding 50 mM Tris pH 8.0 containing 0.1% Tween 20, 0.1 mg/ml lysozyme, 1.0 ug/ml DNAse at room temperature for 30 minutes, followed by centrifugation. Supernatants were then transferred to another 96 well plate and stored at -20 °C until use, when they were thawed at room temperature for 30 minutes. After

shaking briefly, the 2 µl of each culture sample was transferred to a 96-well PCR plate (Bio-Rad, Hercules, CA, USA) containing 7 µl samples buffer (Caliper) followed by heating the plate to 90 °C for 5 minutes on a thermostatically controlled plate heater. The plate was allowed to cool before adding 35 µl water to each sample. The plate was then placed in the instrument along with a protein standard supplied and calibrated by the manufacturer. The instrument functions by mixing the sample with a fluorescent dye that attaches non-covalently to the proteins, followed by electrophoresis through a gel matrix. As the proteins move past a focal point in the chip, the fluorescence signal is recorded and the protein concentration is determined by quantitating the signal relative to the signal generated by a calibrated set of protein standards.

Table 11-7. k_{cat} and K_M and specific activity values for isoprene synthases

Isoprene Synthase	$k_{cat} \pm$ S.D. (s^{-1})	$K_M \pm$ S.D.(mM)	$K_i \pm$ S.D.(mM)	S.A. (nmol/mg/min)
Full Length	0.72±0.09	2.4±0.3	15.7±0.2	420±60
Truncated	1.5±0.2	1.8±0.2	9.8±0.9	800±100
MD09-165	0.8±0.1	2.4±0.5	19±4	440±80
MD09-167	1.1±0.3	2±1	8.7±0.9	610±60

[0382] Parameters were determined by fitting the following uncompetitive substrate inhibition equation to data obtained for rate/[E] vs. [DMAPP]:

$$\frac{rate}{[E]} = \frac{k_{cat} * [S]}{K_M + [S] * \left(1 + \frac{[S]}{K_i}\right)}$$

[0383] The specific activities (S.A.) were calculated for reactions containing 5 mM DMAPP, 50 mM MgCl₂, 100 mM Tris, 100 mM NaCl, and 2.5-4.5 µg isoprene synthase from the supernatant of whole cell lysate. Reactions were performed at 30 °C for 15 minutes in triplicate using independently grown cultures.

Results/Discussion:

[0384] The specific activity of each protein was determined (Figure 40 and Table 11-7). The specific activity of truncated isoprene synthase was approximately 2-fold greater than the specific activity of full length isoprene synthase. The C-terminally His-tagged full length isoprene synthase yielded approximately the same specific activity as the full length isoprene

synthase. The C-terminally His-tagged truncated isoprene synthase gave a specific activity that was less than the non-tagged truncated isoprene synthase, but greater than the specific activity of both full length isoprene synthases.

[0385] The rate of conversion of DMAPP to isoprene was analyzed over a range of DMAPP concentrations in order to determine the k_{cat} and K_M parameters of the enzymes (Figures 41 and 42 and Table 11-7). The enzymes all exhibited rate profiles consistent with uncompetitive substrate inhibition by DMAPP. The K_i for the full length constructs was greater than the K_i for truncated constructs as determined by altering K_i and observing the best fit (R-value) to the data (data not shown). All data has been fit to a K_i of 10 mM in the figures in this document. The K_M s of truncated isoprene synthase decreased relative to the full length isoprene synthases. Therefore, as the substrate concentration decreases the ratio between the isoprene synthase activity of the truncated isoprene synthase to the full length isoprene synthase will increase (Figure 43). The k_{cats} of truncated isoprene synthases increased relative to the full length isoprene synthases. This results in greater isoprene synthase activity for the truncated isoprene synthases than the full length isoprene synthases at all substrate concentrations (Figure 43). The ratio of the isoprene synthase activity of the truncated isoprene synthase to the full length isoprene synthase at varying DMAPP levels was determined (Figure 43).

[0386] Conclusions: “Truncated” isoprene synthases exhibit increased k_{cat} values and decreased K_M values with respect to the catalysis of the conversion of DMAPP to isoprene compared to “full length” isoprene synthases. The specific activity of “truncated” isoprene synthases is increased relative to the specific activity of “full length” isoprene synthases. The most active isoprene synthase was non-tagged truncated isoprene synthase “MEA” variant (in strain MD09-173). The truncated isoprene synthase may increase the isoprene synthase activity in organisms used for production of isoprene compared to the full length isoprene synthase.

Example 12

Isoprene Synthase Enrichment – Relief of DMAPP Toxicity

[0387] This example identifies residue changes within IspS that confer better activity to the enzyme through kinetic improvement, increased expression, increased solubility, or any other means by which DMAPP is more effectively converted to isoprene by Isoprene Synthase. This procedure allows for the relief of cytotoxic intracellular levels of DMAPP by expression of improved IspS variants. In a population of cells expressing a heterogeneous mixture of

IspS variants, the best enzymes should allow for better growth of their host strain, and should be enriched in the mixed population.

I. Strain Construction

[0388] Construction of Strains MCM518-521 and 528-531 in which Lambda promoters drive integrated mKKDyI was as follows. Primers MCM120 and MCM224 (Table 12-1) were used to amplify the resistance cassette from the GeneBridges FRT-gb2-Cm-FRT template using Stratagene Herculanase II Fusion kit according to the manufacturer's protocol. Four 50 μ L reactions were cycled as follows: 95°C, 2 min; (95 °C, 20 sec, 55 °C, 20 sec, 72 °C, 1 min) for 30 cycles; 72 °C, 3 min; and 4°C until cool. The four reactions were pooled and purified on a Qiagen PCR column according to the manufacturer's protocol and eluted with 60 μ L EB at 55 °C.

[0389] Plasmid pRedET-carb (GeneBridges) was electroporated into MCM446. Transformants were recovered by shaking for one hour in SOC (Invitrogen) at 30 °C and then selected on LB containing carbenicillin (50 μ g/ml) (carb50) plates at 30 °C overnight. A carbenicillin resistant colony was frozen as MCM508 (Table 12-2).

[0390] Strain MCM508 was grown from a fresh streak in 5mL LB/carb50 at 30 °C to an OD₆₀₀ of ~0.5. 40mM L-arabinose was added and culture was incubated at 37 °C for 1.5 hrs. Cells were harvested and electroporated with 3 μ L of purified amplicons as previously, and then recovered in 500 μ L SOC at 37 °C for 1.5-3hrs. Transformants were selected on LB/kan10 plates at 37 °C.

[0391] Recombination of the amplicon at the target locus was confirmed by PCR with primers GB-DW and MCM208. The resulting amplicons were sequenced to identify four clones with the sequences below. Carbenicillin-sensitive clones were frozen as strains MCM518-521.

[0392] MCM518-521 were restreaked on LB kan10 plates and grown overnight at 37 °C.

[0393] Strains MCM518-521 were cultured in LB/kan10 at 37 °C and then electrotransformed with plasmid pCP20 (Datsenko KA, Wanner BL. One-step inactivation of chromosomal genes in *Escherichia coli* K-12 using PCR products. Proc Natl Acad Sci U S A. 2000 Jun 6;97(12):6640-5). Cells were recovered in 500 μ L SOC, shaking at 30 °C for 1 hour. Transformants were selected on LB/carb50 plates at 30 °C overnight. The following morning a colony from each transformation was grown at 30 °C in liquid LB/carb50 until visibly turbid. The culture was then shifted to 37 °C for at least 3 hrs. Cells were streaked from this culture onto LB plates and grown overnight at 37 °C.

[0394] The following day colonies were patched to LB, LB/carb50 and LB/kan10. Clones that grew on neither carb50 nor kan10 and were cultured in liquid LB from the patch on LB and frozen as MCM528-531.

DNA Sequences

[0395] These assemblies include the new promoters inserted on the chromosome in strains MCM518-521, as well as the very beginning of the mMVK ORF. Upstream of these assemblies is sequence from the GeneBridges FRT-gb2-Cm-FRT cassette. Downstream is the remainder of the mMVK ORF and then the rest of the lower MVA pathway integron from strain MCM508.

MCM518

```
aaagaccgaccaagcgacgtctgagagctccctggcgaattcggtagcaataaaagagcttt
atcttcatgatctgtgtgttggtttttgtgtgctggcgcggaagttcctattctctagaaagt
ataggaacttcctcgagccctatagtgagtcgtattaaattcatataaaaaacatacagata
accatctgcggtgataaattatctctggcggtgttgacataaataccactggcggtgatact
gagcacatcagcaggacgcactgaccaccatgaaggtgcaaaggaggtaaaaaacatggta
tcctgttctgcgccgggtaagatttacctgttcggtgaacacgccgtagtttatggcgaaac
tgcaattgcgtgtgcggtggaactgcgtaccctgttcgcgcggaactcaatgactctatca
ctattcagagc (SEQ ID NO:56)
```

MCM519

```
aaagaccgaccaagcgacgtctgagagctccctggcgaattcggtagcaataaaagagcttt
atcttcatgatctgtgtgttggtttttgtgtgctggcgcggaagttcctattctctagaaagt
ataggaacttcctcgagccctatagtgagtcgtattaaattcatataaaaaacatacagata
accatctgcggtgataaattatctctggcggtgttgacataaataccactggcggtgatact
gagcacatcagcaggacgcactgaccaccatgaaggtgcaaaggaggtaaaaaacatggta
tcctgttctgcgccgggtaagatttacctgttcggtgaacacgccgtagtttatggcgaaac
tgcaattgcgtgtgcggtggaactgcgtaccctgttcgcgcggaactcaatgactctatca
ctattcagagc (SEQ ID NO:57)
```

MCM520

```
aaagaccgaccaagcgacgtctgagagctccctggcgaattcggtagcaataaaagagcttt
atcttcatgatctgtgtgttggtttttgtgtgctggcgcggaagttcctattctctagaaagt
```

ataggaacttcctcgagccctatagtgagtcgtattaaattcatataaaaaacatacagata
 accatctgcggtgataaattatctctggcggtgttgacctaataaccactggcggtgatact
 gagcacatcagcaggacgcactgaccaccatgaaggtgcaaaggtaaaaaaacatggatatcc
 tgttctgcgccgggtaagatttacctgttcggtgaacacgccgtagtttatggcgaaactgc
 aattgcgtgtgcggtggaactgcgtaccggtgttcgcgcggaactcaatgactctatcacta
 ttcagagc (SEQ ID NO:58)

MCM521 (in strains MCM531 and MD09-171)

aaagaccgaccaagcgacgtctgagagctccctggcgaattcggtagccaataaaagagcttt
 attttcatgatctgtgtgttggtttttgtgtgcggcgcggaagttcctattctctagaaagt
 ataggaacttcctcgagccctatagtgagtcgtattaaattcatataaaaaacatacagata
 accatctgcggtgataaattatctctggcggtgttgacgtaataaccactggcggtgatact
 gagcacatcagcaggacgcactgaccaccatgaaggtgcaaaggaggtaaaaaacatggta
 tcctgttctgcgccgggtaagatttacctgttcggtgaacacgccgtagtttatggcgaaac
 tgcaattgcgtgtgcggtggaactgcgtaccggtgttcgcgcggaactcaatgactctatca
 ctattcagagc (SEQ ID NO:59)

[0396] The neo-PL.2-mKKDyI (from MCM521) was transduced into BL21(λDE3) to generate strain MD09-171. A P1 lysate of MCM521 was made and transduced into BL21(λDE3) according to standard molecular biology techniques (Miller, A Short Course in Bacterial Genetics). Transductants were selected on Kan20 LB medium plates. Positive colonies were further verified by PCR to confirm the presence of PL.2-mKKDyI in the BL21(λDE3) strain. 1 µl of pCP20 plasmid was then transformed into this strain and selected for on LB + Carb50 and incubated at 30 °C. Positive transformants were subsequently streaked on an LB plate and incubated at 37 °C for loss of the pCP20 plasmid. To confirm the loss of the neomycin (kanamycin) resistance marker, colonies that grew at 37 °C were patched onto LB Kan20, LB Carb50, and plain LB plates. The strains with integrated PL.2 mKKDyI without the kanamycin resistance marker that have lost pCP20 should be sensitive to kanamycin and carbenicillin. 4 KanS CarbS were used to check by PCR for the presence of mKKDyI in BL21(λDE3) with the parental BL21(λDE3) strain as a control. Once PCR confirmed the presence of mKKDyI, the resulting strain was transformed with 1 µl of the pLysS plasmid (Invitrogen). The resulting strain, MD09-171, was used for the enrichment experiments described below.

Table 12-1. Primers used for strain construction

MCM120	aaagtagccgaagatgacggtttgtcacatggagttggcaggatgtttgat taaaagcaattaaccctcactaaagggcgg (SEQ ID NO:60)
MCM208	gctctgaatagtgatagagtca (SEQ ID NO:61)
MCM224	taaatcttaccggcgcgacagaacaggataccatgtttttttacctcctttgc accttcatggtggtcagtgcgtcctgctgatgtgctcagtatcac cgccagtggatatttangtcaacaccgccagagataatttatcaccgcagat ggttatctgtatgttttttatatgaatttaatacgactcactatag ggctcg (SEQ ID NO:62)
GB-DW	aaagaccgaccaagcgacgtctga (SEQ ID NO:63)
MCM161	caccatgggtatcctgttctgcg (SEQ ID NO:64)
MCM162	ttaatctactttcagaccttgc (SEQ ID NO:65)
MCM143	aggaggtggtctcaaatagactgccgacaacaatagta (SEQ ID NO:66)
MCM144	aggaggtggtctcagcgctctgcagttatagcattctatgaatttgctg (SEQ ID NO:67)

Table 12-2. Strains

Strain	Description	Parent
MCM508	BL21 gi1.6-mKKDyI + predet.-carb	MCM446
MCM518	BL21 neo-PL.6-mKKDyI, clone10	MCM508
MCM519	BL21 neo-PL.0-mKKDyI, clone11	MCM508
MCM520	BL21 neo-PL.0-mKKDyI (bad RBS in front of mMVK), clone13	MCM508
MCM521	BL21 neo-PL.2-mKKDyI, clone15	MCM508
MCM528	BL21 PL.6-mKKDyI, loopedout	MCM518
MCM529	BL21 PL.0-mKKDyI, loopedout	MCM519
MCM530	BL21 PL.0-mKKDyI (bad RBS in front of mMVK), loopedout	MCM520
MCM531	BL21 PL.2-mKKDyI, loopedout	MCM521
MD09-171	BL21 (λ DE3) PL.2-mKKDyI, loopedout + pLysS	MCM521

II. Growth Inhibition of MCM531 by Mevalonic Acid

[0397] An overnight culture of MCM531 (see strain description) was back-diluted to an OD₆₀₀ of 0.05 (this corresponds to an OD₆₀₀ of approx 0.005 in a 96-well plate reader) (SpectraMax M2, Molecular Devices). The diluted culture was then aliquotted into separate wells in a 96-well deep-well plate into standard TM3 medium (13.6 g K₂PO₄, 13.6 g KH₂PO₄, 2.0 g MgSO₄*7H₂O) supplemented with 1% glucose and 0.8 g/L Biospringer yeast extract (1% Yeast extract final))with 0, 1, 5, 10, 15 or 20 mM MVA added. Figure 44 shows the growth curve of MCM531 in the various concentrations of MVA. Each MVA concentration was assayed in quadruplicate, error bars were negligible. Figure 44 shows that

MCM531 was severely inhibited for growth at concentrations of 5 mM MVA and higher.

III. Mutagenesis of IspS and Selection/Enrichment Assay

[0398] To generate a randomly mutagenized IspS open reading frame, the GeneMorph II EZ Clone domain mutagenesis kit (Stratagene) was used according to the manufacturer's recommended protocol. Specific primers to amplify the template (Pdu39 (pET24a-P. alba (MEA))) are described below (Table 12-3, pET24 Megaprime Forward and Reverse). To achieve the desired mutation frequency, the protocol outlined in the GeneMorph II kit was followed. To generate 2 to 3 residue changes per molecule, approximately 150 ng of starting template DNA was used for the initial PCR reaction. More or less template was used to generate IspS enzymes with either fewer or more residue changes, respectively. The resulting mutant "megaprimers" were then used to amplify the rest of the plasmid according to the manufacturer's recommended protocol.

[0399] The final PCR product from the GeneMorph II kit was treated with *DpnI* according to the manufacturer's recommended protocol. Prior to transformation into *E. coli*, it was necessary to desalt the PCR reaction by microdialysis. Typically, approximately 20 μ l of the PCR reaction was subjected to microdialysis and used for electroporation into strain MD09-171 (Table 12-2) by standard molecular biology procedures. After electroporation, cells were recovered for 2 hours at 30 °C, and then plated onto LB medium Kan50 Cm35. The entire transformation volume was plated to recover all mutations generated by the mutagenesis procedure.

[0400] For enrichment, all transformants were scraped and pooled together. An aliquot from each pool was frozen for storage in the -80. Prior to the enrichment, strains (control or experimental pools) were inoculated directly into LB with Kan50 and grown for a few hours, to minimize the loss of pool heterogeneity. After this recovery period, cells were diluted into liquid TM3 medium (described above) with Kan50, 5 mM MVA, and 200 μ M IPTG. (The exact dilution was determined empirically for each pool/source of medium/control reaction). Cultures were then placed in the shaking incubator at 34°C until an OD₆₀₀ of no more than 5. At this point, plasmids were purified from 1 ml of the "enriched" culture via miniprep protocol (Qiagen). This plasmid preparation was then transformed into electrocompetent MD09-171 cells as described above. The transformed cells were recovered for 2 hours in LB medium without antibiotics, and then subjected to a subsequent round of enrichment by dilution into TM3 with Kan50, 5 mM MVA and 200 μ M IPTG as described above. This culture was placed into the shaking incubator at 34°C until it reached an OD₆₀₀ of no more

than 5, as described above. Plasmids were then purified, retransformed and subjected to further rounds of “enrichment.” The enrichment process continued for 5 or 6 rounds of selection, plasmid purification, and retransformation. The process continued until the culture was homogeneous, *i.e.* contained only one variant of IspS by sequencing analysis.

[0401] After the last round of enrichment, the plasmid pool was transformed into chemically competent *E. coli* Top10 cells (Invitrogen) per the manufacturer’s recommended protocol, recovered, plated onto LB medium with Kan50, and sent for complete sequencing (Quintara Biosciences) for comparison to the wild type sequence of *P. alba* IspS. Primers used for sequencing are described below.

Table 12-3. Primers used for IspS mutagenesis and sequencing

pET24 Megaprime Forward	gtttaactttaagaaggagatatatacat
pET24 Megaprime Reverse	gagctcgaattcggatcctta
alba sequencing reverse	ctcgtacaggctcaggatag
alba sequencing reverse2	ttacgtcccaacgctcaact
EWL1000	gcactgtcttccgtctgctgc
QB1493	cttcggcaacgcatggaaat

IV. Individual Residue Changes Identified by Enrichment/DMAPP Toxicity Relief:

[0402] Plasmids isolated from selection/enrichment were fully sequenced. The following residue changes were confirmed by sequencing (Quintara Biosciences). The residue numbering corresponds to the *P. alba* “Full Length” sequence (in *P. alba* pET24a; SEQ ID NO:120), where the starting methionine is amino acid number 1. Identified residue changes included: V10M, F12S, T15A, E18G, V58I, V58F, L70Q, L70V, L70T, T71P, V79L, E89D, G94A, S119F, F120L, G127R, E175V, T212I, S257A, R262G, A266G, F280L, N297K, F305L, L319M, E323K, A328T, D342E, A359T, K366N, E368D, L374M, S396T, V418S, K438N, H440R, T442I, T442A, I449V, A469S, K500R, K505Q, G507S, S509N, F511Y, and N532K.

[0403] Combinations of Residue Changes (in a single IspS enzyme) Identified by Enrichment/DMAPP Toxicity Relief: G127R/F511Y, L70Q/G94A/R262G/F305L, F12S/T15A/E18G/N297K, S396T/T442I, V10M/E323K, F120L/A266G, K438N/K500R, V79L/S509N, E175V/S257A/E368D/A469S, T71P/L374M, F280L/H440R, E89D/H440R, V58F/A328T/N532K, S119F/D342E/I449V, and K366N/G507S.

Example 13

Construction of Site Saturation Libraries (SSLs) and Biochemical Analysis of the L70R Variant

[0404] This example includes an examination of possible amino acid substitutions at sites identified by the selection/enrichment procedure described above and other sites of potential interest (active site, conserved between Poplar species) for their effect on solubility, expression, and activity of IspS.

I. Strain Construction

[0405] Residues identified by the selection/enrichment procedure (L70, G94, R262, F305) described above were chosen for analysis. In addition, residues that are putatively involved in substrate binding (F303, V3065, F385, S412, Q416, F450), and residues that are different between the various Poplar species (*e.g.* V418, T442) were chosen for further analysis (numbering corresponds to the “full length” amino acid sequence of *P. alba* IspS). To generate a randomized pool of amino acid substitutions (the Site Saturation Library, SSL), pDu39 (see description above) was subjected to QuickChange (Stratagene) mutagenesis with the QC primers indicated below (Table 13-1), according to the manufacturer’s recommended protocol. The PCR Reaction was prepared as follows:

1 µl pDu39, 5 µl 10X PfuUltra HF buffer, 1 µl dNTP’s, 1 µl (50 µM) primer-For (*e.g.* QC L69 F), 1 µl (50uM) primer-Rev (*e.g.* QC L69 R), 2 µl DMSO, 39 µl diH₂O, and 1 µl PfuUltra HF Polymerase (Stratagene). The PCR cycling parameters for QuickChange were as follows: 95°C 1 min, 95°C 30 sec., 55°C 1 min., 68°C 7.3 min. for one cycle followed by 95°C 30 sec., 55°C 1 min., 68°C 7.3 min for 17 more cycles. The temperature was then reduced to 4°C. Incorporation of the bases NNK at the codon for the residues described above allows for the insertion of codons that represent all 20 possible amino acids at the given site. 1 µl of the resulting pools of mutagenized PCR products was *DpnI* treated (described above), and transformed into chemically competent *E. coli* Top10 cells (Invitrogen) according to the manufacturer’s recommended protocol. Entire transformation reactions were recovered in 1 ml LB medium without antibiotics at 37 °C for 1 hour and plated onto LB Kan50. The next day, all transformants were scraped off of the LB plates, mixed thoroughly, and plasmids were purified via miniprep (Qiagen). Pools of plasmids were then transformed into chemically competent BL21(λDE3)pLysS cells (Invitrogen) according to the manufacturer’s protocol. The transformation reactions were recovered in 1 ml LB medium at 37 °C for 1 hour and then plated onto LB Kan50 Cm35 at dilutions sufficient to generate separation of positive colonies. After overnight incubation at 37 °C,

individual colonies were inoculated into individual wells in a 96-well deep-well microtiter plate (VWR) containing 500 μ l of liquid LB Kan50 Cm35 each. In eight wells (typically column 12, A through H) strain MD09-173 (see above) or BL21(λ DE3)pLysS with P. alba pET24a (Full Length) was inoculated as a control for the DMAPP assay. The microtiter plates were then sealed with a semi-permeable membrane (Breathe-Easier, Diversified Biotech), and incubated overnight at 30 °C in a shaking incubator (Vertiga). The next day, 100 μ l samples from each well within a 96-well plate were mixed with 50 μ l of 50% glycerol in a new 200 μ l 96-well plate, and frozen at -80 °C until further analysis. This plate was then used for the DMAPP assay described below.

Table 13-1. Primers used for mutagenesis

QC L70 F	gaaaaagcagaattttnnkaccctgctggaactg (SEQ ID NO:68)
QC L70 R	cagttccagcagggtmnnaaattctgctttttc (SEQ ID NO:69)
QC G94 F	gagtctgatatccgtnnkgcgctggatcgcttc (SEQ ID NO:70)
QC G94 R	gaagcgatccagcgcmmnacggatatcagactc (SEQ ID NO:71)
QC R262 F	tcccgttggtggcgtnnkgtgggtctggcgacc (SEQ ID NO:72)
QC R262 R	ggtcgccagacccacmnnacgccaccaacggga (SEQ ID NO:73)
QC F303 F	tccgtcgcaaaaatggnktctttcgtaaccatt (SEQ ID NO:74)
QC F303 R	aatggttacgaaagamnncaatttttgcgacgga (SEQ ID NO:75)
QC F305 F	gcaaaaatgttttctnnkgtaaccattatcgac (SEQ ID NO:76)
QC F305 R	gtcgataatggttacmnnagaaaacatttttgc (SEQ ID NO:77)
QC V306 F	aaaatgttttctttcnnkaccattatcgacgat (SEQ ID NO:78)
QC V306 R	atcgtcgataatggtmnnagaaagaaaacatttt (SEQ ID NO:79)
QC F385 F	gacctgtgcaacgctnnkctgcaagaagccaag (SEQ ID NO:80)
QC F385 R	cttggcttcttgcagmnnagcgttgacaggtc (SEQ ID NO:81)
QC S412 F	gcatggaaatcctctnnkggcccgcgtgcaactg (SEQ ID NO:82)
QC S412 R	cagttgcagcgggcmnnagaggatttccatgc (SEQ ID NO:83)
QC Q416 F	tcttctggcccgcgtggnkctggtgttcgcttac (SEQ ID NO:84)
QC Q416 R	gtaagcgaacaccagmnnacgcgggccagaaga (SEQ ID NO:85)
QC V418 F	ggcccgcgtgcaactggnkttcgttacttcgct (SEQ ID NO:86)
QC V418 R	agcgaagtaagcgaamnnacagttgcagcgggcc (SEQ ID NO:87)
QC T442 F	caaaaataccatgacnnkatctctcgtccttcc (SEQ ID NO:88)
QC T442 R	ggaaggacgagagatmnngtcatggtatttttg (SEQ ID NO:89)
QC F450 F	cgtccttcccatatcnnkcgctctgtgcaatgac (SEQ ID NO:90)
QC F450 R	gtcattgcacagacgmnnagatatgggaaggacg (SEQ ID NO:91)

II. Generation of a “Winner” Plate for Secondary Assay and Identification of L70R as a Variant with Increased Specific Activity

[0406] Variants that displayed increased specific isoprene production when compared to wild type were chosen for further analysis. Figure 45 shows a typical data set of an SSL plate for an individual residue, in this case L70. From this particular plate, the samples in wells C3 (27), D3 (39), and E3 (51) were chosen for further analysis. Other variants at different residues (listed above) that showed increased isoprene productivity, when subjected to DMAPP analysis, were isolated from their original SSL plates stored at -80 °C (described above), and re-arrayed onto a new “winner” plate for secondary screening. Two wells containing MD09-173 were included as controls. All variants were sequenced (Quintara Biosciences) and subjected to the DMAPP assay as described above. See Table 13-2 for sequencing results. For the DMAPP assay, a single growth of each variant was assayed, and therefore a single lysate generated, but in quadruplicate to generate statistically significant data. Samples were assayed at the OD₆₀₀ indicated in Table 13-2. Protein analysis was performed on all lysates using the Western Breeze Western blot kit (Invitrogen) followed by fluorescence detection on a Storm860 (see below)

III. DMAPP Assay – Growth, Lysis and Isoprene Measurement

[0407] A patch plate was prepared from a glycerol stock plate using a VP-Scientific Replication Tool patch LB Agar CM35/Kan50 large patch plate from overnight glycerol stock Plate. Cultures were incubated at 30°C overnight (20 to 24 hr). Plates were stored at 4°C for up to a week.

[0408] An overnight growth plate was prepared from 500 mL of LB CM35/Kan50 media. 300 µL/well of LB CM35/Kan50 media were dispensed into deep 96 well plate. Using the V&S Replication Tool, the patch inoculum was transferred to a deep 96 well plate. Media was inoculated by dipping the tool then shaking the pin within the well. The overnight growth plate was sealed with a Breathe-Easier Sealing Membrane. The plate was incubated at 30 °C overnight in a Vertiga Shaking Incubator at 800 rpm for 16 to 18 hours.

[0409] A deep 96 well day growth plate was prepared by dispensing 588 µL/well of LB CM35/Kan50 media. The overnight growth plate was removed from the incubator and cultures were diluted 50-fold. 12 µL of overnight culture was transferred to a day growth plate containing 588 µL/well of supplemented LB media. The overnight growth plate was sealed with a new Breathe-Easier Sealing Membrane and was incubated at 34 °C and 800 rpm for 2.25 hr in the Vertiga Shaking Incubator.

[0410] To induce the expression of IspS, thawed 12 mM IPTG was poured into 50 mL or 100 mL sterile reservoir and dispensed 20 µL/well into each 600 µL/well culture. Overnight

growth plate was resealed with Breathe-Easier Sealing Membrane and incubated at 34°C and 800 rpm for 4 hours in the Vertiga Shaking Incubator.

[0411] To harvest cell, 200 µL of induced culture was transferred to 450 µL Nunc storage plate. The plate was centrifuged at 3300 rpm for 20 min at 4 °C in a low speed benchtop centrifuge. 180 µL supernatant was removed with a pipettor and discarded. The plate was sealed with an aluminum foil membrane, covered with a plastic plate lid and stored frozen at -80 °C.

[0412] The OD₆₀₀ of the plates were read. 150 µL 1x PBS was dispensed into a 96 well Costar Read Plate (#9017). 50 µL of culture sample was then transferred to the read plate. The OD₆₀₀ reading was then taken with a Spectramax Plate Reader.

[0413] Lysis: The harvest plate was defrosted in a room temperature water bath for 4 min and then incubated in Thermomixer at 25°C at 1200 rpm for 1 min. Lysis buffer was dispensed at 80 µL/well to 20 µL/well of harvest cells. 1.25X Lysis Working Stock buffer was prepared from 6.25 ml 1M Tris pH 8, 625 µl 10% Tween 20, 312.5 µl 0.2 M PMSF, 462.5 µl 10 mg/ml DNase I (Sigma), 1.25 ml 1 M MgCl₂, 132.5 µl 25000 U/µl Lysozyme (Epicentre Technologies) and 40.968 ml dH₂O. Plates were incubated on a Thermomixer at 25°C and 1200 rpm for 30 min. 1x Lysis Buffer stock was prepared by diluting 19 ml of 1.25x Lysis Working stock with 4.75 ml dH₂O.

[0414] DMAPP Working Stock was prepared as follows.

Dimethylallyl Pyrophosphate (triammonium salt) 25 mg (Cayman Chemicals, Cat No. 63180)

final	total vol	diluent	stock vol	stock	dilution	Reagent	Units	Plate(s)
30	840	840	solid	25000	1	DMAPP	mg/mL	1
3	2500	2250	250	30	10	DMAPP	mg/mL	1

[0415] The Diluent was 0.1 M Potassium Phosphate. 200 uL/well was dispensed for transfer to sample wells and was stored on ice.

[0416] DMAPP Reaction: 1x Lysis Buffer was dispensed at 65 µL/well. 15 µL/well of lysate was transferred to the respective sample wells in a 96 Deep well Zinsser Glass Block. DMAPP reagent was dispensed at 20 µL/well. The glass block was sealed with an aluminum foil membrane and incubated at 25 °C at 450 rpm for 45 min. The reaction was stopped by transferring the Glass Block to a 70 °C water bath and incubating for 6 min. GC Analysis was performed as previously described.

IV. Western blot of Isoprene Synthase with Fluorescence labeled Secondary Antibody.

[0417] Samples were prepared and run on NativePAGE™ Novex® Bis-Tris Gels (Invitrogen) according to the manufacturer's protocol. After completion of the run the gels were immediately transferred to Nitrocellulose membranes using the XCell II™ Blot Module (Invitrogen) according to the manufacturer's recommended protocol. After transfer, the membrane was placed in 15 ml of the appropriate Blocking Solution (Ultra filtered Water 31.5 ml, Blocker/Diluent (Part A) 9 ml, Blocker/Diluent (Part B) 4.5 ml) in the covered, plastic dish provided in the kit and incubated for 30 minutes on a rotary shaker set at 1 revolution/sec. The Blocking Solution was decanted and the membrane rinsed 2 times with 20 ml of water for 5 minutes. The membrane was incubated with 15 ml of Primary Antibody (Ab) Solution (24 µl of primary Ab in 15 ml Blocking solution) for 1 hour, followed by washing 3 times 5 minutes with 20 ml of 1x Antibody Wash Solution. The membrane was then incubated in 15 ml of Secondary Antibody Solution (15 µl secondary Ab (Alexa Fluor 488 goat anti-rabbit IgG (H+L, Invitrogen)) in 15 of blocking solution) for 30 minutes. The membrane was washed 3 times at 5 minutes with 20 ml of Antibody Wash, and rinsed two times 2 minutes with 20 ml of water. The membrane was dried between paper towels and stored at room temperature for further detection. The fluorescent bands were detected and quantified using the Storm 860 Molecular Imager (GMI, Inc).

V. Results:

[0418] Table 13-2 shows all of the relevant data for each variant assayed: sequencing results, residue change, average isoprene production, protein concentration, and average specific activity (of all 4 replicates). Figure 46 shows the graphical representation of the data shown in Table 13-2. Specific activity was calculated by multiplying the isoprene produced (µg/l) by 0.0414 and then dividing by protein concentration (mg/ml). This conversion factor (0.0414) accounts for the total headspace volume in a sealed 2 ml GC vial (1.9 ml), the lysate volume (15 µl), the duration of the DMAPP assay (45 min), and the molecular weight of isoprene. Thus, specific activity values are given in nmol isoprene/mg Isps/min.

[0419] The data in Table 13-2 and the graph in Figure 46 show that of all variants analyzed, all three L70R variants displayed higher specific activity than wild type. To analyze the L70R variants further, the specific activity values for all 3 isolates (4 replicates of each) were averaged and compared to the controls (2 isolates, 4 replicates each). Therefore, there were 12 measurements for the L70R variant, and 8 for wild type. These data are shown

below in Table 13-3. When corrected for protein, the L70R variants displayed a 25% increase in activity over MD09-173 (the MEA truncation). Figure 47 shows the average specific activity for all L70R variants compared to MEA. Error bars show one standard deviation. The two data sets for L70R and the MEA control were subjected to a Student's T-Test for statistical analysis, which yielded a P-value of 6.0011×10^{-5} .

Table 13-2. Sequencing Results, Isoprene Production, Protein concentration, and Specific Activity for all residues in the Winner Plate

Residue	WT codon	Mutant codon	Amino acid change	OD600	Avg. Isoprene produced (4 replicates)	Std. Dev. (isoprene)	Protein concentration (mg/ml)	Avg. Specific Activity (4 replicates)	Std. Dev. (Specific Activity)
L70	CTG	TGG	W	2.6488	663.00	18.64	0.139	198.020	5.566
L70	CTG	CGG	R	3.4944	1235.72	48.82	0.127	401.451	15.860
L70	CTG	CGT	R	3.32416	944.77	19.15	0.121	322.523	6.539
L70	CTG	CGT	R	3.37232	947.11	17.74	0.106	371.375	6.957
L70	CTG	TGG	W	2.83584	774.49	32.60	0.129	248.925	10.477
G94	GGT	GAG	E	3.304	1070.41	20.52	0.144	307.731	5.898
G94	GGT	GGG	G	3.40816	982.64	27.37	0.133	306.891	8.549
R262	CGT	CGT	R	3.01616	861.08	29.13	0.116	307.227	10.395
R262	CGT	CGG	R	3.31744	965.44	19.25	0.136	294.311	5.868
F305	TTC	TTT	F	2.90864	788.46	15.15	0.119	274.297	5.270
F305	TTC	CTG	L	2.7776	626.70	22.91	0.141	183.823	6.719
WT1				3.71392	991.65	35.78	0.144	285.216	10.292
F305	TTC	TTC	F	2.604	694.03	35.93	0.118	244.488	12.658
V306	GTA	GTA	V	3.33648	912.56	38.76	nd	nd	nd
V306	GTA	GTA	V	3.17632	892.97	56.10	0.137	269.415	16.927
F385	TTC	TTC	F	2.89184	721.24	29.36	0.134	223.380	9.094
Q416	CAG	CAG	Q	3.00272	953.07	57.38	0.149	264.718	15.938

V418	GT G	ATG	M	2.1414 4	767.25	27.94	0.164	194.09 9	7.067
V418	GT G	ACG	T	3.2468 8	973.57	75.08	0.192	209.80 1	16.17 9
V418	GT G	ACG	T	3.1326 4	936.98	68.50	0.171	226.19 7	16.53 6
T442	AC C	GTA	V	3.0755 2	951.64	38.06	0.154	255.49 6	10.21 7
F450	TTC	TTC	F	2.9332 8	911.44	50.88	0.126	299.01 8	16.69 1
V418	GT G	ACG	T	3.0889 6	1057.71	46.28	0.183	238.88 7	10.45 4
WT2				3.3756 8	960.25	58.16	0.131	304.34 6	18.43 4

[0420] Note that sequence of all plasmids is identical to Pdu39 (see above) with the exception of the indicated codon. The L70R variants are highlighted in gray.

Table 13-3. Average Specific Activity of all L70R variants relative to the MEA control.

VariantAverage	Specific Activity	Standard Deviation
L70R	365.116	35.31977
WT (the MEA control)	294.7809	17.19228

See Figure 47 for bar graph.

Example 14

Truncations of *P. alba*, *P. tremuloides*, *P. trichocharpa*, and Kudzu

Isoprene Synthases

[0421] This example describes the generation of a series of truncations in the IspS enzymes of *P. alba*, *P. tremuloides*, *P. trichocharpa*, and Kudzu and to determine their effect on activity.

I. Strain Construction

[0422] All isoprene synthase genes were codon optimized for *E. coli*, synthesized, and cloned into pET24a by DNA2.0 (Menlo Park, CA). All truncated constructs were generated using the QuickChange Site-Directed Mutagenesis kit (Stratagene) using the previously described templates *P. alba* pET24a (for plasmids pDu47-3 through -7, Figures 48, 49, 51-60), *P. tremuloides* pET24a (plasmid pDu48, Figures 49C, 61 and 62), *P. trichocharpa* pET24a (pDu49, Figures 50A, 63, 64), or pET24d-Kudzu (pDu50 and 50-4, Figures 50B, 50C, 65-68) for PCR amplification. Approximately 50 ng of template DNA was used for

amplification (with an Eppendorf Mastercycler Gradient PCR Machine) of the PCR product with the Forward (For) and Reverse (Rev) primer pairs that correspond to each relevant truncation (QC Trunc -3 F and QC Trunc -3 R, for example, see Table 14-1). PCR reactions mixtures were as follows: 1 μ l P. alba pET24a (or other template), 5 μ l 10X PfuUltra HF buffer, 1 μ l dNTP's (10 mM), 1 μ l (50uM) primer-For, 1 μ l (50 μ M) primer-Rev, 1.5 μ l DMSO, 39.5 μ l diH₂O and 1 μ l PfuUltra HF Polymerase. PCR cycle parameters were as follows: (95 °C 1 min., 95 °C 1 min., 55 °C 1 min., 68 °C 7.30 min.) for 18 cycles then 4°C until cool using an Eppendorf Mastercycler Gradient Machine. The PCR products were treated with 1-2 μ l of *DpnI* (Roche) for 1-3 hour at 37 °C. 5 μ l of the *DpnI* treated products were visualized on a 0.8% E-gel (Invitrogen). 1 μ l of each product was transformed into chemically competent *E. coli* Top10 cells (Invitrogen) (according to the manufacturer's protocol). Transformants were selected for on LB medium containing kanamycin at a concentration of 50 μ g/ml (Kan50), and incubated overnight at 37 °C. Five colonies of each transformation were selected and grown to stationary phase in 3 ml liquid LB Kan50. Plasmids were purified using a QIAprep Spin miniprep kit (Qiagen) according to the manufacturer's recommended protocol. Purified plasmids were sequenced (by Quintara Biosciences) with T7 Forward and Reverse primers, compared to the parental sequence, and confirmed for their respective truncation. The resulting plasmids (pDu47-3 through pDu50-4, see Table 14-2) were transformed into chemically competent *E. coli* BL21(DE3) pLysS (Invitrogen) according to the manufacturer's recommended protocol. Table 14-3 describes the strains used for expression of truncated IspS enzymes.

[0423] After overnight incubation at 37 °C, individual colonies were inoculated into individual wells in a 96-well deep-well microtiter plate (VWR) containing 500 μ l of liquid LB Kan50 CM35 each. Microtiter plates were then sealed with a semi-permeable membrane (Breathe-Easier, Diversified Biotech), and incubated overnight at 30 °C in a shaking incubator (Vertiga). The next day, 100 μ l samples from each well within a 96-well plate were mixed with 50 μ l of 50% glycerol in a new 200 μ l 96-well plate, and frozen at -80 °C until further analysis. This plate was then used for the DMAPP assay as described in Example 13. Table 14-4 shows the average specific productivity of all samples, and Figure 69 shows the graphical representation of the same data.

[0424] DMAPP activity and protein quantitation was determined as described in Example 13.

[0425] Specific activity was calculated by multiplying the isoprene produced (μ g/l) by 0.00776 and then dividing by protein concentration (mg/ml). This conversion factor

(0.00776) accounts for the total headspace volume in a sealed 2 ml GC vial (1.9 ml), the lysate volume (80 μ l), the duration of the DMAPP assay (45 min), and the molecular weight of isoprene. Thus, specific activity values are given in nmol isoprene/mg IspS/min.

Table 14-1. Primers

QC Trunc -3 F	gaaggagatatacatatgaccgaagctcgtcgt (SEQ ID NO:92)
QC Trunc -3 R	acgacgagcttcggtcatatgtatatctccttc (SEQ ID NO:93)
QC Trunc -4 F	gaaggagatatacatatggaaaccgaagctcgt (SEQ ID NO:94)
QC Trunc -4 R	acgagcttcggtttccatatgtatatctccttc (SEQ ID NO:95)
QC Trunc -5 F	gaaggagatatacatatgactgaaaccgaagct (SEQ ID NO:96)
QC Trunc -5 R	agcttcggtttcagtcatatgtatatctccttc (SEQ ID NO:97)
QC Trunc -6 F	gaaggagatatacatatggaaactgaaaccgaa (SEQ ID NO:98)
QC Trunc -6 R	ttcggtttcagtttccatatgtatatctccttc (SEQ ID NO:99)
QC Trunc -7 F	gaaggagatatacatatgaccgaaactgaaacc (SEQ ID NO:100)
QC Trunc -7 R	ggtttcagtttcggtcatatgtatatctccttc (SEQ ID NO:101)
QC Kudzu MEA F	agaaggagatataccatggaagctcgtcgttcgcaaac (SEQ ID NO:102)
QC Kudzu MEA R	gtttgcggaacgacgagcttccatggtatatctccttct (SEQ ID NO:103)
QC Kudzu -4 F	agaaggagatataccatggagcataattcccgt (SEQ ID NO:104)
QC Kudzu -4 R	acgggaattatgctccatggtatatctccttct (SEQ ID NO:105)
QC Trem/Trich -2 F	gaaggagatatacatatggaaacgcgctcgttct (SEQ ID NO:106)
QC Trem/Trich -2 R	agaacgacgcgtttccatatgtatatctccttc (SEQ ID NO:107)

Table 14-2. Plasmids

pDu47-3	Mtg pET24a-P.alba TRC (-3)
pDu47-4	Mtg pET24a-P.alba TRC (-4)
pDu47-5	Mtg pET24a-P.alba TRC (-5)
pDu47-6	Mtg pET24a-P.alba TRC (-6)
pDu47-7	Mtg pET24a-P.alba TRC (-7)
pDu48	Mtg pET24a-P.tremu TRC (MET)
pDu49	Mtg pET24a-P.tricho TRC (MET)
pDu50	Mtg pET24d-Kudzu TRC (MEA)
pDu50-4	Mtg pET24d-Kudzu TRC (-4)

Table 14-3. Strains

MD09-197-3	BL21(DE3)pLysS, pDu47-3
MD09-197-4	BL21(DE3)pLysS, pDu47-4
MD09-197-5	BL21(DE3)pLysS, pDu47-5
MD09-197-6	BL21(DE3)pLysS, pDu47-6
MD09-197-7	BL21(DE3)pLysS, pDu47-7
MD09-198	BL21(DE3)pLysS, pDu48
MD09-199	BL21(DE3)pLysS, pDu49
MD09-200	BL21(DE3)pLysS, pDu50
MD09-200-4	BL21(DE3)pLysS, pDu50-4

Table 14-4. Specific Productivity of Variants listed above

Sample	Isoprene Produced (µg/l)				OD ₆₀₀		OD Norm
	Sample	Average	Std Dev	% CV	Raw	xDF xCF	
IspS N terminal Truncation Variants	MD09-197-3	3381	154	5	0.3060	3.4	986.6
	MD09-197-4	3000	206	7	0.3160	3.5	847.7
	MD09-197-5	2932	266	9	0.3200	3.6	818.2
	MD09-197-6	2450	217	9	0.3040	3.4	719.4
	MD09-197-7	2285	397	17	0.3020	3.4	675.5
	MD09-198	1916	106	6	0.3330	3.7	513.6
	MD09-199	2031	108	5	0.2140	2.4	847.6
	MD09-200	141	8	6	0.3700	4.1	34.0
	MD09-200-4	1829	197	11	0.2760	3.1	591.5
	MD09-173	2414	201	8	0.3400	3.8	633.9
	MD09-176	2826	354	13	0.3260	3.7	773.9
	BL21 DE3 pLysS + P alba pET24a	2175	117	5	0.2990	3.3	649.6

[0426]

[0427] A second experiment was conducted with the strains outlined in Table 14-5.

Control was BL21 DE3 pLysS with P. alba pET24a (full length P. alba IspS).

[0428]

[0429] Table 14-5. Strains

MD09-197-3	BL21(DE3)pLysS, pDu47-3
MD09-197-4	BL21(DE3)pLysS, pDu47-4
MD09-197-5	BL21(DE3)pLysS, pDu47-5
MD09-197-6	BL21(DE3)pLysS, pDu47-6
MD09-197-7	BL21(DE3)pLysS, pDu47-7
MD09-198	BL21(DE3)pLysS, pDu48
MD09-199	BL21(DE3)pLysS, pDu49
MD09-173	BL21(DE3)pLysS, pET24a-P.alba (MEA) Untagged (pDu39)
MD09-174	BL21(DE3)pLysS, pET24a-P.alba (MNV) Untagged (pDu40)
MD09-175	BL21(DE3)pLysS, pET24a-P.alba (MSV) Untagged (pDu41)
MD09-176	BL21(DE3)pLysS, pET24a-P.alba (MTE) Untagged (pDu42)
MD09-177	BL21(DE3)pLysS, pET24a-P.alba (MVS) Untagged (pDu43)
MD09-197-3	BL21(DE3)pLysS, pDu47-3
MD09-197-4	BL21(DE3)pLysS, pDu47-4

[0430]

[0431] Results

[0432] All truncations of *P. alba* IspS and two from *P. tremuloides* and *P. trichocharpa* were assayed in parallel to compare their relevant specific activities via DMAPP assay and quantitative Western blot. At least two samples per variant were assayed for isoprene

production and amount of IspS in mg/ml. Variant MD09-174 produced little isoprene and expressed little protein, yet displayed high specific activity. High specific activities were also displayed by MD09-173, MD09-176, and MD09-197-3 (see Table 14-6 and Figure 70). The highest levels of protein (μg) in 3 μg total protein were displayed by MD09-176 and MD09-197-3, indicating that these variants are more effectively expressed in the *E. coli* BL21 DE3 host strain.

[0433]

[0434] Table 14-6. Specific Activity of truncations.

Strain	Average Specific Activity	Standard Deviation	Average ug IspS/ 3 μg total protein
P. alba FL	240.16239	31.0423851	0.198
MD09-175	331.8755329	8.958408319	0.188
MD09-177	340.1959506	39.72104203	0.150
MD09-176	363.516921	3.376026129	0.202
MD09-174	452.7792122	27.71567075	0.018
MD09-197-7	279.2042431	23.82331163	0.158
MD09-197-6	309.8357305	7.903564316	0.164
MD09-197-5	293.22592	20.97161876	0.165
MD09-197-4	321.3926574	56.13760028	0.186
MD09-197-3	333.0604155	45.92710529	0.207
MD09-173	368.4597405	37.80631246	0.159
MD09-198	297.6476631	56.81405985	0.154
MD09-199	256.7342861	8.239697653	0.216

[0435]

Example 15

Constructs for Three-Dimensional Structure Determination

I. Construction of pMAL-C4X Kudzu

[0436] A synthetic gene, coding for isoprene synthase (IspS) of the Kudzu vine (*Pueraria lobata*) and codon-optimized for *E. coli*, was purchased from DNA2.0 (Menlo Park, CA) and provided as plasmid p9795 (Figure 71 and 72). The insert was removed by digestion with *Bsp*LU11I/*Pst*I, gel-purified, and religated into *Nco*I/*Pst*I-digested pTrcHis2B (Invitrogen, Carlsbad, California). The resulting plasmid was named pTrcKudzu (Figure 73: map of pTrcKudzu). The stop codon in the insert is before the *Pst*I site, which results in a construct in which the His-Tag is not attached to the IspS protein.

[0437] A PCR reaction was performed to amplify the *E. coli* codon-optimized Kudzu gene using plasmid pTrcKudzu as the DNA template, primers EL-959 and EL-960, 10mM dNTP

(Roche, Indianapolis, IN), and Pfu Ultra II Fusion DNA polymerase (Stratagene, La Jolla, CA) according to manufacturer's protocol. PCR conditions were as follows: 95 °C for 2 min (first cycle only), 95 °C for 25 sec, 60 °C for 25 sec, 72 °C for 30 sec, repeat for 28 cycles, with final extension at 72 °C for 1 min. The PCR product was then purified using the QIAquick PCR Purification Kit (Qiagen Inc, Valencia, CA).

[0438] The Kudzu PCR product (1 µg) was digested using *Eco*RI and *Hind*III restriction endonucleases (Roche) according to manufacturer's protocol. The digest was incubated 37 °C for 30 minutes to minimize digestion of the internal *Eco*RI site that is present in the Kudzu gene. The digested PCR fragment was then purified using the QIAquick PCR Purification Kit. The vector pMAL-C4X (0.5 µg) (New England Biolabs, Ipswich, MA; Figures 75 and 76) was digested using *Eco*RI and *Hind*III restriction endonucleases (Roche) according to manufacturer's protocol. The digested vector was then gel purified using the QIAquick Gel Extraction Kit (Qiagen Inc). A DNA ligation reaction was performed using T4 DNA ligase (New England Biolabs) with a 5:1 ratio of digested Kudzu PCR product to digested pMAL-C4X vector according to manufacturer's protocol. An aliquot of the ligation reaction was then transformed into TOP10 chemically competent cells (Invitrogen Corp). Transformants were selected on LA + 50 µg/µl carbenicillin plates.

[0439] Screening of transformants containing the Kudzu gene was performed by picking colonies and performing PCR with primers EL-957 and EL-966 using PuReTaq Ready-To-Go PCR beads (GE Healthcare, Piscataway, NJ) according to manufacturer's protocol. PCR conditions were as follows: 95 °C for 2 min (first cycle only), 95 °C for 30 sec, 50 °C for 30 sec, 72 °C for 40 sec, repeat for 28 cycles, with final extension at 72 °C for 1 min. PCR products were analyzed on a 2% E-gel (Invitrogen Corp) looking for a 600 bp fragment. Colonies containing the correct sized PCR product insert were submitted for DNA sequencing using primers EL-950, EL-951, EL-953, and EL-957. DNA sequencing confirmed the construction of plasmid pMAL-C4X Kudzu (Figures 77 - 79).

Table 15-1 Primer sequences

Primer name	Primer sequence
EL-950	CGGTGAACTGAAAGGTGACGTCC (SEQ ID NO:108)
EL-951	GGACGTTAACGCTATTAACACCCTG (SEQ ID NO:109)
EL-953	CACATCGTCGATCAGCTCCAGC (SEQ ID NO:110)
EL-957	GGTCGTCAGACTGTGCGATGAAGCC (SEQ ID NO:111)
EL-959	GCTTATGAATTCTGTGCGACCTCTTCTCAATTTACTCAG (SEQ ID NO:112)
EL-960	GCTTATAAGCTTAGACATACATCAGCTGGTTAATCGGG (SEQ ID NO:113)

EL-966	CTCCTCCAGCAGGTTCTCACC (SEQ ID NO:114)

[0440] Plasmid pMAL-C4X Kudzu was transformed into OneShot BL21(λDE3) chemically competent cells (Invitrogen Corp). Expression strain transformants were selected on LA + 50 mg/ml carbenicillin plates.

II. IspS Variants for Crystal Structure Trials

[0441] This example describes methods to generate affinity tagged isoprene synthase (IspS) enzymes for expression, purification and crystallization.

Strain Construction

[0442] For constructs in the pET200D-TOPO vector (Invitrogen), PCR products of the IspS enzymes from *P. alba*, *P. tremuloides*, and *P. trichocarpha* were gel extracted and purified (Qiagen), using 0.8% E-gel (Invitrogen), according to the manufacturer's recommended protocol. PCR reactions for pET200 constructs are as follows: Reaction mixture was 1 μl (Templates)-pET24a-P.alba, 5 μl 10X PfuUltraII Fusion buffer, 1 μl dNTP's (10 mM), 1 μl primer (50uM) primer F-(MCM219 or 218), 1 μl primer (50uM) primer R-(MCM182), 41 μl diH₂O and 1 μl of PfuUltra II Fusion DNA Polymerase from Stratagene; Cycle Parameter were 95°C 1 min., 95°C 1 min, 55°C 20 sec., 72°C 27 sec. for 29 cycles followed by 72°C for 3 min and then 4° C until cool, using an Eppendorf Mastercycler. Similar reactions were performed for *P. tremuloides*, *P. trichocarpha*, and Kudzu. 3 μl of purified product was then ligated to the pET200D/TOPO vector (Invitrogen), according to the manufacturer's protocol. The reaction was incubated for 5 minutes at room temperature, and the 6 μl topoisomerase mixture was then transformed into E. coli Top10 chemically competent cells (Invitrogen) according to the manufacturer's protocol. Transformants were selected for on LB Kan50, and incubated at 37 °C overnight. Five colonies per construct were chosen and screened using PuReTaq Ready-To-Go PCR Beads (Amersham) using the T7 Forward and MCM182 primers (Table 15-2). Clones harboring inserts of the correct size were further verified by plasmid miniprep using the QIAprep Spin Miniprep kit (Qiagen) followed by sequencing using the T7 Forward and T7 Reverse primers (Quintara Biosciences). One fully sequenced construct for each IspS variant (see below for details and sequence/Figures 79-90), was chosen for further study. 1 μl of each plasmid was transformed into BL21(λDE3) pLysS (Invitrogen) according to the manufacturer's protocol. Transformants were selected for on LB medium with Kan50 + Cm35 and incubated at 37°C overnight. The resulting strains were used for expression and purification of various IspS enzymes for crystallography studies.

[0443] Construction of N-terminally 6His-tagged IspS plasmids, strains and purification is described in Example 11.

Table 15-2. Primers

MCM219	caccatgcgtttagcgtgtcca (SEQ ID NO: 114)
MCM182	gggcccggttaactttaactagactctgcagtttagcgttcaaacggcagaa (SEQ ID NO:115)
MCM218	caccatgcgtcgttctgcgaactac (SEQ ID NO:116)

Table 15-3. Plasmids

P. alba pET24a	pET24a with “full length” IspS from P. alba
P. trichocharpa pET24a	pET24a with “full length” IspS from P. trichocharpa
P. tremuloides pET24a	pET24a with “full length” IspS from P. tremuloides
MBP-Kudzu	
pDu27	P. albaFL-pET200/Top 10
pDu30	P. albaTRC-pET200/Top10
pDu31	P. tremTRC-pET200/Top 10
pDu32	P. trichTRC-pET200/Top10
MD09-161	pET24a-P.alba FL C-Term (+) TEV, His tag/MCM331
MD09-163	pET24a-P.alba TRC (MEA) C-Term (+) TEV, His tag/MCM331

Table 15-4. Strains

MBP-Kudzu	
MD08-99	BL21 DE3 pLys + pDu27
MD08-100	BL21 DE3 pLys + pDu30
MD08-102	BL21 DE3 pLys + pDu31
MD08-104	BL21 DE3 pLys + pDu32
MD09-165	BL21(DE3)pLysS, pET24a-P.alba FL C-Term (+) TEV, His tag
MD09-167	BL21(DE3) pLysS, pET24a-P.alba TRC (MEA) C-Term (+) TEV, His tag

III. Digestion of TEV (Tobacco Etch Virus) or EK (Enterokinase)-tagged enzymes

TEV Cleavage (IspS from strains MD09-165 and MD09-167)

[0444] Strains MD09-165 and MD09-167 are described in Example 11. For digestion, enzymes were purified through a Ni charged sepharose (GE Healthcare) and desalted into 50 mM HEPES, 50 mM NaCl pH 7.4 buffer containing 1 mM DTT. Digestion was performed with TurboTEV Protease from Eton Bioscience Inc. One unit of TurboTEV per 10 µg of purified protein was used. The digest was performed at 4°C overnight. Samples were passed through another Ni column equilibrated in the Ni buffer to remove uncleaved enzyme, tag, TurboTEV protease (that is also tagged) and impurities. The Ni column pass though and washes were analyzed using SDS-PAGE gel (NUPAGE, Invitrogen) and DMAPP activity

assays. Samples containing pure enzyme were pooled and desalted into 50 mM NaCl pH 7.4 buffer containing 1 mM DTT and stored at -80°C.

EK Cleavage (IspS from strains MD08-102 and MD08-104)

[0445] For digestion enzymes were purified through a Ni charged sepharose (GE Healthcare) and desalted into 50 mM HEPES, 50 mM NaCl pH 7.4 buffer containing 1 mM DTT.

[0446] Digestion was performed with EKMax (E180-02) (Invitrogen) using 1 unit of EKMax per 20 µg of purified protein at 4 °C overnight. Samples were passed over EK Away resin (Invitrogen) to remove excess enterokinase. Samples were batched onto Ni charged sepharose resin (equilibrated in the Ni wash buffer) and incubated for 30 min at 4 °C, with occasional inverting. This removed uncleaved enzyme, tag, and impurities. The Ni column pass through and washes were analyzed using SDS-PAGE gel (4-12% NUPAGE, Invitrogen) and DMAPP activity assays. Samples containing pure enzyme were pooled and desalted into 50 mM HEPES, 50 mM NaCl pH 7.4 buffer containing 1 mM DTT and stored at -80°C.

IV. Purification of MBP-IspS

[0447] Construction of pMAL-C4X Kudzu for the expression of MBP-Kudzu isoprene synthase is described above. MBP-Kudzu isoprene synthase production from *E. coli* grown in batch culture at the 15-L scale.

Medium Recipe (per liter fermentation medium)

[0448] K₂HPO₄ 7.5 g, MgSO₄ * 7H₂O 2 g, citric acid monohydrate 2 g, ferric ammonium citrate 0.3 g, yeast extract 0.5 g, 1000X Modified Trace Metal Solution 1 ml. All of the components were added together and dissolved in diH₂O. This solution was autoclaved. The pH was adjusted to 7.0 with ammonium hydroxide (30%) and q.s. to volume. Glucose 10 g, thiamine * HCl 0.1 g, and antibiotics were added after sterilization and pH adjustment.

1000X Modified Trace Metal Solution

[0449] Citric Acids * H₂O 40 g, MnSO₄ * H₂O 30 g, NaCl 10 g, FeSO₄ * 7H₂O 1 g, CoCl₂ * 6H₂O 1 g, ZnSO₄ * 7H₂O 1 g, CuSO₄ * 5H₂O 100 mg, H₃BO₃ 100 mg, NaMoO₄ * 2H₂O 100 mg. Each component is dissolved one at a time in diH₂O, pH to 3.0 with HCl/NaOH, then q.s. to volume and filter sterilized with 0.22 micron filter.

[0450] Fermentation was performed in a 15-L bioreactor with BL21 (DE3) *E. coli* cells containing the pMAL-C4X plasmid expressing a maltose binding protein (MBP)-Kudzu

isoprene synthase fusion molecule. This experiment was carried out to produce isoprene synthase at the desired fermentation pH 7.0 and temperature 30 °C. A frozen vial of the *E. coli* strain was thawed and inoculated into tryptone-yeast extract medium. After the inoculum grew to OD 1.0, measured at 550 nm (OD₅₅₀), 120 mL was used to inoculate a 15-L bioreactor bringing the initial volume to 9-L.

[0451] Expression of the desired molecule was achieved by adding isopropyl-beta-D-1-thiogalactopyranoside (IPTG). The IPTG concentration was brought to 1 mM when the OD₅₅₀ reached a value of 10. Cells containing the desired product were harvested 3 hrs after IPTG addition.

MBP-IspS Purification

[0452] The broth was centrifuged for 15 min at 10000 x g. The pellet was collected and frozen at -80 °C until further purification. Cells were resuspended in MBP-Bind Buffer (5% glycerol, 20 mM Tris pH 7.4, 200 mM NaCl, 2 mM DTT, 1mg/ml lysozyme) and passed through the french press three times at 20000 psi. The lysate was then ultracentrifuged at 100000 x g for 1 hour to yield a relatively clear solution. The supernatant was pipetted from the top of the tube without disturbing the gelatinous material on the bottom of the centrifuge tube. Gel filtration was performed on the supernatant using a Superdex-200 26/60 column (GE healthcare). The column was developed using MBP-Bind buffer at a flow rate of 3 mL/min at 23 °C. Fractions were tested for DMAPP activity as described below. Active fractions were pooled and loaded onto 25 mL amylose resin (New England Biolabs). The column was washed with 10 column volumes MBP-Bind buffer and the protein was then eluted with 2 column volumes of MBP-Bind buffer containing 10 mM maltose to yield >90% pure MBP-IspS.

V. DMAPP Assay

[0453] The following reaction mixture was used for the DMAPP assay: 25 µL lysate mixture, 5 µL MgCl₂ (1 M), 5 µL DMAPP (100mM), and 65 µL 100 mM Tris pH 8, 100 mM NaCl for a total volume of 100 µL. The reaction is performed at 30 °C for 15 minutes in a gas tight 1.8 mL GC tube. Reactions are terminated by the addition of 100 µL 500 mM EDTA (pH 8). The amount of isoprene produced was measured by GC/MS as described above.

Example 16

Three-dimensional structure of IspS

[0454] Seven constructs of plant isoprene synthase (IspS) were prepared to generate crystals suitable for x-ray diffraction. These were: a construct containing N-terminal histidine-tagged maltose binding protein and kudzu IspS (MBP-kudzu), full-length *P. alba* IspS with N-terminal histidine-tag (MD08-99), *P. alba* IspS with the first nineteen N-terminal residues removed (MD08-100), this construct also had the N-terminal his-tag removed after purification. Full-length, untagged *P. alba* IspS (strain RM11608-2). A truncated *P. alba* IspS construct featuring two additional residues before the twin-arginine motif was generated (MD09-167). *P. tricharpa* IspS was generated, which contains both an N-terminal his-tag and N-terminal truncation (MD08-104), and another construct composed of IspS from *P. tremuloides* was generated with an N-terminal his-tag and N-terminal truncation (MD08-102). Construction of strains expressing various isoprene synthases are described above.

[0455] Each construct was purified and a concentrated protein solution was then prepared for surveying possible crystallization conditions. Each construct was purified independently and surveyed as described below. All in-house crystallization screens were set up using the hanging drop vapor diffusion method. At a minimum, each construct was surveyed using the following commercial screens: the Crystal Screen from Hampton Research (Aliso Viejo, CA) and the JCSG+ Suite from Qiagen (Valencia, CA).

[0456] Purified MBP-kudzu using was set up using the following commercial screens: the Crystal Screen from Hampton Research and the JCSG+ Suite from. Additionally, purified MBP-kudzu was sent to the Hauptman-Woodward Institute (Buffalo, NY) for high-throughput screening, where no fewer than 1536 conditions were surveyed. The purified MBP-kudzu fusion precipitated out of solution in the majority of conditions, and no protein crystals were observed.

[0457] The next construct used for crystallization screening was MD08-99 (full-length *P. alba* IspS with N-terminal histidine-tag). MD08-99 was purified and the histidine-tag was removed. The same three initial crystallization screens were performed as for MBP-kudzu. The Hampton Research Crystal Screen and Qiagen JCSG+ Suite were each performed at multiple protein concentrations. Small needle-like crystals were observed in some Hampton Research Crystal Screen conditions. Further attempts to improve the crystals involved co-crystallization with the IspS inhibitor sodium ibandronate (Sigma-Aldrich, St Louis, MO). Taken together, an additional 288 crystallization conditions were attempted with variations of pH, concentration, and crystallization reagents. The nine best crystals were then prepared for data collection and tested in-house on a Rigaku RU200 rotating anode generator and R-Axis IV++, and they either did not diffract x-rays or were salt crystals.

[0458] The first nineteen N-terminal residues of *P. alba* IspS were removed to produce construct MD08-100. This construct had the N-terminal histidine-tag removed after purification. In house crystallization screens were performed using the Hampton Research Crystal Screen and Qiagen JCSG+ Suite, each with multiple protein concentrations. Initial crystal hits included hexagonal plates that diffracted to 16Å resolution, and small rods that diffracted to 5Å resolution using the in-house x-ray generator. In an attempt to improve the crystals, MD08-100 was co-crystallized with either sodium ibandronate or sodium pyrophosphate (Sigma-Aldrich, St Louis, MO), both of which are inhibitors of the IspS activity. Neither inhibitor resulted in improved crystals or improved diffraction. An additional 168 crystallization conditions were attempted with variations of pH, concentration, and crystallization reagents. The twenty-one most promising MD08-100 crystals were screened for diffraction, with the best resolution obtained being 5Å.

[0459] Full-length, untagged *P. alba* IspS (strain RM11608-2) from a fermentation run was purified. An initial screen was set up using the Hampton Research Crystal Screen, and crystals were observed in four different conditions. All four crystals were tested for diffraction in-house, with three being salt crystals and one not diffracting.

[0460] A truncated *P. alba* IspS construct featuring two additional residues before the twin-arginine motif was generated (MD09-167). This construct contains a C-terminal histidine-tag, and crystallization experiments were set up with the tag either cleaved or not cleaved, at varying protein concentrations, and with or without sodium pyrophosphate. Initial crystallization screens were done as per MBP-kudzu. Crystals from this construct were observed in numerous conditions; optimization included 528 variations of pH, precipitating agents, concentrations, and inhibitors. From the optimization experiments, fifteen different MD09-167 crystals were screened in-house for diffraction. In an effort to improve the resolution, various crystal freezing conditions were tested, with the effect of improving the diffraction limits from 10Å to 6.5Å.

[0461] A new construct containing *P. tricharpa* IspS was generated, which contains both an N-terminal histidine-tag and an N-terminal truncation (MD08-104). Purified MD08-104 with cleaved histidine-tag was surveyed using the Hampton Research Crystal Screen and the Qiagen JCSG+ suite. This construct generated heavier precipitate than the *P. alba* IspS constructs. Very small needles were observed, with none of the crystals being suitable for diffraction.

[0462] Another construct composed of IspS from *P. tremuloides* was generated with an N-terminal histidine-tag and an N-terminal truncation (MD08-102). Purified MD08-102 with

and without cleaved histidine-tag was set up using the Hampton Research Crystal Screen and the Qiagen JCSG+ Suite at varying protein concentrations. Rod and plate-like crystals were observed in some conditions and an additional 120 experiments were performed to improve the crystals by varying pH, concentration, and crystallization reagents. From the optimization experiments, ten crystals were tested in-house, with the initial best diffraction reaching 5Å. Upon further modification of the freezing conditions of the crystals, a crystal was found that diffracted to 3.3Å from the non-cleaved histidine-tagged protein. This crystal was grown by mixing 2 µL of protein (10 mg/ml, with 30mM MgCl₂) with 2 µL of precipitant solution [10% (wt/vol) polyethylene glycol 8000, 0.1 M HEPES, pH 7.5, 8% ethylene glycol] and equilibrated against 500 µL of precipitant. A cluster of rod-shaped crystals appeared after three weeks. The crystals belong to the tetragonal space group P43212, and have unit cell dimensions $a = 154.2$, $b = 154.2$, $c = 142.7$.

[0463] In-house x-ray diffraction data were collected under a nitrogen stream at 100 K using a Rigaku RU200 generator and R-Axis IV++ detector. Before flash-freezing the crystal in liquid nitrogen, it was cryoprotected by swiping it through a solution containing 10% (wt/vol) polyethylene glycol 8000, 0.1 M HEPES, pH 7.5, and 25% ethylene glycol. Data were integrated using Mosflm (Leslie, A. (1998) *J. of Appl. Crystallography* 30, 1036-1040) and scaled using SCALA (Collaborative Computational Project, N. (1994) *Acta Crystallographica Section D* 50, 760-763). The data were then phased by molecular replacement using MrBUMP (Keegan, R. M., and Winn, M. D. (2007) *Acta Crystallographica Section D* 63, 447-457; Vagin, A., and Teplyakov, A. (1997) *J. of Appl. Crystallography* 30, 1022-1025), with a monomer of limonene synthase (Protein Data Bank ID 2ONH)(Berman, H., *et al.* (2007) *Nucl. Acids Res.* 35, D301-303) as the starting model. The crystal contains one dimer in the asymmetric unit with a solvent content of 66%.

[0464] A 3.05Å data set from the same crystal was then collected using beamline 11-1 of the Stanford Synchrotron Radiation Laboratory. These data were also processed using Mosflm and SCALA. Data collection and refinement statistics are given in Table 16-1.

[0465] Refinement with Refmac5 (Collaborative Computational Project, N. (1994) *Acta Crystallographica Section D* 50, 760-763) was used with iterative manual rebuilding steps using the visualization program Coot (Emsley, P., and Cowtan, K. (2004) *Acta Crystallographica Section D* 60, 2126-2132). During refinement, the geometry of the protein was checked using Molprobity (Davis, I. W., *et al.* (2007) *Nucl. Acids Res.*, gkm216).

[0466] The fold of poplar IspS is similar to bornyl diphosphate synthase (Whittington, D. A., *et al.* (2002) *Proc. Natl Acad. Sci. USA* 99, 15375-15380), limonene synthase (Hyatt, D.

C., et al. (2007) Proc. Natl Acad. Sci. USA 104, 5360-5365), and tobacco 5-*epi*-aristolochene synthase (Starks, C. M., *et al.* (1997) Science 277, 1815-1820). The structure consists of two helical domains, a C-terminal domain containing the active site and N-terminal domain (Figures 90 and 91). Coordinates are provided in Table 16-7.

Table 16-1. Data Collection and Refinement Statistics

Data Collection	
Space Group	P4 ₃ 2 ₁ 2
Cell dimensions	
a, b, c (Å)	154.2, 154.2, 142.7
α, β, γ, (°)	90.0, 90.0, 90.0
Resolution (Å)	37.8 – 3.05
R _{merge}	16.4 (72.9) ^a
<I/sI>	10.3 (2.6)
Completeness (%)	99.8 (100)
Redundancy	7.3 (7.4)
Refinement	
Resolution (Å)	37.8 – 3.05
No. measured reflections	248741
No. Unique reflections	34201
R _{work}	21.1
R _{free}	27.1
rmsd bonds, (Å)	0.011
rmsd angles, (°)	1.28
No. of Atoms	
Protein, ions ^b	8331
Water	18

Flexible Loops

[0467] The unique and unexpected discovery coming from the determination of the three dimensional structure of isoprene synthase is that several crucial loops forming the active site are flexible. The discovery can be immediately seen when the known structure of other terpene synthases are compared with the structure of isoprene synthase (Figure 94). Overall, the structures are highly conserved in the conformation of secondary structure and connectivity loops. (In this example of the *Poplar tremuloides* IspS from construct P.tremTRC-pET200, the numbering convention is such that the first number of the complete sequence containing the tag is -35, with the first residue of IspS being 1.) However, three segments, forming a considerable portion of the substrate binding pocket, notably the truncated N-terminus, along with two loops comprised of residues 438-453 (Loop I) and

residues 512-527 (Loop II) are seen to diverge (Figures 95 to 97). This has been attributed to the absence of substrate complexed with the enzyme in our structure determination.

[0468] In comparing the enzyme with BdpS, for example, we find that the loops corresponding to residues 498-513 and 573-587 are composed of the same number of residues and have a homologous, but not identical amino sequence in these regions. We expect that the related terpene synthases will be found to display similar flexibility in the segments as these structure become more thoroughly studied. The residues in terpene synthases corresponding to these variable loop regions are enumerated in Table 16-2.

Table 16-2. Residues corresponding to variable loops in terpene synthases

	Poplar IspS	LS	BdpS	TEAS
N-term I	Met 1	Met 57	Ile 54	Val 14
	Arg 2	Arg 58	Arg 55	Arg 15
	Arg 3	Arg 59	Arg 56	Pro 16
	Ser 4	Ser 60	Ser 57	Val 17
	Ala 5	Gly 61	Gly 58	Ala 18
	Asn 6	Asn 62	Asn 59	Asp 19
	Tyr 7	Tyr 63	Tyr 60	Phe 20
	Glu 8	Asn 64	Gln 61	Ser 21
	Pro 9	Pro 65	Pro 62	Pro 22
	Asn 10	Ser 66	Ala 63	Ser 23
	Ser 11	Arg 67	Leu 64	Leu 24
	Trp 12	Trp 68	Trp 65	Trp 25
	Asp 13	Asp 69	Asp 66	Gly 26
	Tyr 14	Val 70	Ser 67	Asp 27
	Asp 15	Asn 71	Asn 68	Gln 28
	Tyr 16	Phe 72	Tyr 69	Phe 29
N-term II	Leu 17	Ile 73	Ile 70	Leu 30
	Leu 18	Gln 74	Gln 71	Ser 31
	Ser 19	Ser 75	Ser 72	Phe 32
	Ser 20	Leu 76	Leu 73	Ser 34
	Asp 21	Leu 77	Asn 74	Ile 35
	Thr 22	Ser 78	Thr 75	Asp 36
	Asp 23	Asp 79	Pro 76	Asn 37
	Glu 24	Tyr 80	Tyr 77	Gln 38
	Ser 25	Lys 81	Thr 78	Val 39
	Ile 26	Glu 82	Glu 79	Ala 40
	Glu 27	Asp 83	Glu 80	Glu 41
	Val 28	Lys 84	Arg 81	Lys 42

Loop I	Leu 438	Leu 498	Leu 498	Thr 446
	Ala 439	Gly 499	Gly 499	Ala 447
	Ser 440	Thr 500	Thr 500	The 448
	Ala 441	Ser 501	Ser 502	Tyr 449
	Ser 442	Val 502	Tyr 503	Glu 450
	Ala 443	Glu 503	Phe 504	Val 451
	Glu 444	Glu 504	Glu 505	Glu 452
	Ile 445	Val 505	Leu 506	Lys 453
	Ala 446	Ser 506	Ala 507	Ser 454
	Arg 447	Arg 507	Arg 508	Arg 455
	Gly 448	Gly 508	Gly 509	Gly 456
	Glu 449	Asp 509	Asp 510	Gln 457
	Thr 450	Val 510	Val 511	Ile 458
	Ala 451	Pro 511	Pro 512	Ala 459
	Asn 452	Lys 512	Lys 513	Thr 460
	Ser 453	Ser 513	Thr 514	Gly 461
Loop II	Tyr 512	Tyr 576	Tyr 573	Tyr 520
				Ile 521
	His 513	His 577	Leu 574	His 522
	Asn 514	Asn 578	His 575	Asn 523
	Gly 515	Gly 579	Gly 576	Leu 524
	Asp 516	Asp 580	Asp 577	Asp 525
	Ala 517	Gly 581	Gly 578	Gly 526
	His 518	His 582	Phe 579	Tyr 527
	Thr 519	Gly 583	Gly 580	The 528
	Ser 520	Thr 584	Val 581	His 529
	Pro 521	Gln 585	Gln 582	Pro 530
	Asp 522	His 585	His 583	Glu 531
	Glu 523	Pro 586	Ser 584	Lys 532
	Leu 524	Ile 587	Lys 585	Val 533
	Thr 525	Ile 588	Thr 586	Lue 534
	Arg 526	His 589	Tyr 587	Lys 535

[0469] This important finding can be exploited for the engineering of improved isoprene synthase in a straightforward manner. It would be desirable to exploit the flexibility to enhance enzyme performance by making substitutions in the amino acids forming these segments to facilitate the transitions the enzyme must undergo in the steps of binding substrate and allowing rearrangement of substrate in different kinetic steps that are postulated

to occur during enzymatic de-phosphorylation and for electron transfer to convert DMAPP to isoprene.

[0470] The structure provides the new insight that these loops can be present in at least two conformations: the “open” form in the absence of substrate, as we have see in the uncomplexed structure of the isoprene synthase, and a “closed,” or active form when the substrate is bound. It would therefore also be beneficial to modify residues coming in contact with the loops in the active form as described in Table 16-3.

Table 16-3. Residues coming within 5 angstroms of flexible elements

	<i>P trem</i> IspS	1N1B	2ONG	5EAS ^c
N-term neighbors	L17, L18, S19, S20, S239, R243, F253, A254, R255, D256, R257, I259, E260, D293, Y295, D296, V297, Y298, G299, T300, E303, Y325, L374, Y375, elements of loop I, elements of loop II, V529, L530, T534	70I, 71Q, 72S, 73L, 298S, 302S, 312F, 313V, 314R, 315D, 316R, 318V, 319E, 352D, 354Y, 355D, 356V, 357Y, 358G, 359T, 362E, 384Y, 433Y, 434H, elements of loop I, elements of loop II, 589I, 590A, 594F	I73, Q74, S75, L76, F299, R303, F313, A314, R315, A316, R317, V319, E320, D353, Y355, D356, V357, Y358, G359, T360, E363, Y385, F434, Y435, elements of loop I, elements of loop II, M590, T591, F595	L30, S31, F32, S33, S248, K252, Y262, A263, R264, D265, R266, V268, E269, D302, F304, D305, A306, Y307, G308, T309, E312, Y334, F383, I384, elements of loop I, elements of loop II, I538, I538, V543
Loop I neighbors	Elements of N-term, D293, Y295, V297, E370, A371, W373, L374, S378, T379, P380, F382, Y385, F386, R433, L434, C435, N436, D437, V454, S455, C456, Y457, M458, T469, V472, I476, Y512, elements of loop II	Elements of N-term, 352D, 354Y, 356V, 429E, 430A, 432W, 433Y, 437Y, 438T, 439P, 441L, 444Y, 445L, 493R, 494L, 495P, 496D, 497D, 514I, 515Q, 516C, 517Y, 518M, 529V, 532V, 536I, 572Y, elements of loop II	Elements of N-term, D353, Y355, V357, E430, A431, W433, F434, H438, K439, P440, L442, Y445, L446, R493, L494, A495, D496, D497, L514, Q515, C516, Y517, M518, R529, V532, I536, Y573, elements of loop II	Elements of N-term, D304, F304, A306, E379, S380, W382, F383, Y387, T388, P389, V391, Y394, L395, R441, V442, I443, D444, D445, I462, E463, C464, C465, M466, M477, F480, A484, Y520, elements of loop II
Loop II neighbors	Elements of N-terminus, E187, L188, R255, R257, F270, E271, Q273, Y274, F285, V288, A439, S440, S442, S508, H509, C510, T511, Y512, R528, V529, L530, S531, V532	Elements of N-terminus, 246D, 247L, 314R, 316R, 329E, 330S, 332F, 333W, 344I, 348I, 499G, 500T, 503Y, 568A, 569Q, 570F, 571I, 572Y, 588H, 589I, 590A, 591G, 592L	Elements of N-terminus, D247, I248, R315, R317, E330, P331, Q333, H334, N345, I347, G499, T500, V502, A569, Q570, L571, M572, Y573, Q589, M590, T591, R592, T593	Elements of N-terminus, E195, Q196, R264, R266, F279, E280, Q282, Y283, I294, I297, A447, T448, E450, V516, E517, V518, T519, Y520, H537, I538, I539, N540, L541

Selection of Sites for Improvement of Plant Isoprene Synthase

[0471] The isoprene synthases of plants were expected to be homologous to the terpene synthases. The three-dimensional structures of three homologous terpene synthases have been determined: *Salvia officinalis* bornyl diphosphate synthase (BdpS; pdb entry 1N1B), *Mentha spicata* limonene synthase (LS; pdb entry 2ONG), and tobacco 5-epi-aristolochene synthase (TEAS; pdb entry 5EAS). These enzymes share only 33% homology but their tertiary structure is conserved. Sequence identity is shown in Table 16-4, and structural homology between the structures is shown in Table 16-5. In addition, the structures of intermediate complexes with all three related enzymes have shown that not only tertiary

folding, but also detailed interactions in the active sites of these enzymes are highly conserved.

Table 16-4. Percent Identity of Terpene Cyclases.

	<i>P trem</i> IspS	1N1B	2ONG	5EAS ^e
<i>P alba</i> IspS ^a	98.6	40.7	41.3	33.2
<i>P trem</i> IspS ^b		41.0	41.4	33.2
1N1B ^c			51.4	33.8
2ONG ^d				33.3

^a Polar alba isoprene synthase

^b Polar tremuloides isoprene synthase

^c bornyl diphosphate synthase

^d limonene synthase

^e 5-epi-aristolochene synthase

Table 16-5. Structural Alignment of Terpene Synthases

	1N1B ^b	2ONG ^c	5EAS ^d
<i>P trem</i> IspS ^a	1.40 (465) ^e	1.29 (468)	1.62 (458)
1N1B		1.27 (520)	1.97 (476)
2ONG			1.83 (477)

^a *Polar tremuloides* isoprene synthase

^b bornyl diphosphate synthase

^c limonene synthase

^d 5-epi-aristolochene synthase

^e Root mean square deviation in Å for Ca atoms, with the number of aligned residues in parenthesis

[0472] In this example of the *Poplar tremuloides* IspS from construct P.tremTRC-pET200, the numbering convention is such that the first number of the complete sequence containing the tag is -35, with the first residue of IspS being 1.

[0473] A comparison of the active site from the structure of BdpS and the structure of poplar IspS indicates that the active site involved in metal ion binding and phosphate recognition is conserved. In particular, Arg 255, Asp 292, Asp 296, Glu 370, Arg 433 and

Asn 436 of poplar IspS were observed to overlap equivalent residues in BdpS. The positioning of an intermediate of the BdpS was also compared with the poplar IspS structure. Based on this it was possible to identify the analogous binding region and the approach direction that dimethylallyl pyrophosphate would require in order to bind and react with the poplar IspS enzyme.

[0474] Based on the structure of poplar IspS, sites in the poplar IspS were identified as candidates for mutagenesis to produce variant IspS enzymes with improved performance. Briefly, sites were selected in the IspS that might alter the interaction of the metal binding, diphosphate recognition, DMAPP chain binding and/or the approach to the active site.

I. Diphosphate/ Metal Binding Sites

[0475] The side chains of amino acid residues in the poplar IspS that are found in proximity to the metal and diphosphate binding side chains were identified. These residues include Asp 293, Tyr 385, Ser 392, and Asp 437. Engineering of these sites may result in increased enzyme activity.

II. Substrate Access Loops

[0476] The substrate access loops of poplar IspS are in regions that deviate from the BdpS structure. In the BdpS structure the residues create a cover over the active site. It is likely that upon substrate binding the structure of poplar IspS will form a similar structure. As such the residues in these loops, including residues 440-453 and 512-524, may be in a position to alter the activity of the poplar IspS. In the poplar IspS enzyme, residues 440-453 have the sequence SASAEIARGETANS and residues 512-526 have the sequence YHNGDAHTSPDEL.

III. Isoprenyl Binding Site

[0477] The complex of BdpS and the product of the reaction, bornyl diphosphate (PDB entry 1N24), was used to identify residues in the poplar IspS structure that with protein engineering may be used modulate substrate specificity and/or reaction rate (altered on and off rates of substrate and product). These residues include Ser 261, Trp 264, Phe 285, Thr 289, Ser 393, Ser 394, Phe 432, and Try 512.

Table 16-6. Candidate mutagenesis sites.

	Poplar IspS
DPP/Metal Binding Sites	Asp 293
	Tyr 385
	Ser 392
	Asp 437
Substrate Access Loop I	Ser 440
	Ala 441
	Ser 442
	Ala 443
	Glu 444
	Ile 445
	Ala 446
	Arg 447
	Gly 448
	Glu 449
	Thr 450
	Ala 451
	Asn 452
	Ser 453
Substrate Access Loop II	Tyr 512
	His 513
	Asn 514
	Gly 515
	Asp 516
	Ala 517
	His 518
	Thr 519
	Ser 520
	Pro 521
	Asp 522
	Glu 523
	Leu 524
Isoprenyl Binding Site	Ser 261
	Trp 264
	Phe 285
	Thr 289
	Ser 393
	Ser 394
	Phe 432

	Tyr 512

Table 16-7. Coordinates of *P. tremuloides* IspS

```

HEADER ---- XX-XXX-XX xxxx
COMPND ---
REMARK 3
REMARK 3 REFINEMENT.
REMARK 3 PROGRAM : REFMAC 5.5.0088
REMARK 3 AUTHORS : MURSHUDOV,VAGIN,DODSON
REMARK 3
REMARK 3 REFINEMENT TARGET : MAXIMUM LIKELIHOOD
REMARK 3
REMARK 3 DATA USED IN REFINEMENT.
REMARK 3 RESOLUTION RANGE HIGH (ANGSTROMS) : 3.05
REMARK 3 RESOLUTION RANGE LOW (ANGSTROMS) : 110.17
REMARK 3 DATA CUTOFF (SIGMA(F)) : NONE
REMARK 3 COMPLETENESS FOR RANGE (%) : 99.67
REMARK 3 NUMBER OF REFLECTIONS : 32446
REMARK 3
REMARK 3 FIT TO DATA USED IN REFINEMENT.
REMARK 3 CROSS-VALIDATION METHOD : THROUGHOUT
REMARK 3 FREE R VALUE TEST SET SELECTION : RANDOM
REMARK 3 R VALUE (WORKING + TEST SET) : .21396
REMARK 3 R VALUE (WORKING SET) : .21092
REMARK 3 FREE R VALUE : .27112
REMARK 3 FREE R VALUE TEST SET SIZE (%) : 5.1
REMARK 3 FREE R VALUE TEST SET COUNT : 1727
REMARK 3
REMARK 3 FIT IN THE HIGHEST RESOLUTION BIN.
REMARK 3 TOTAL NUMBER OF BINS USED : 20
REMARK 3 BIN RESOLUTION RANGE HIGH : 3.050
REMARK 3 BIN RESOLUTION RANGE LOW : 3.129
REMARK 3 REFLECTION IN BIN (WORKING SET) : 2359
REMARK 3 BIN COMPLETENESS (WORKING+TEST) (%) : 100.00
REMARK 3 BIN R VALUE (WORKING SET) : .288
REMARK 3 BIN FREE R VALUE SET COUNT : 127
REMARK 3 BIN FREE R VALUE : .352
REMARK 3
REMARK 3 NUMBER OF NON-HYDROGEN ATOMS USED IN REFINEMENT.
REMARK 3 ALL ATOMS : 8349
REMARK 3
REMARK 3 B VALUES.
REMARK 3 FROM WILSON PLOT (A**2) : NULL
REMARK 3 MEAN B VALUE (OVERALL, A**2) : 24.592
REMARK 3 OVERALL ANISOTROPIC B VALUE.
REMARK 3 B11 (A**2) : .41
REMARK 3 B22 (A**2) : .41
REMARK 3 B33 (A**2) : -.81
REMARK 3 B12 (A**2) : .00
REMARK 3 B13 (A**2) : .00
REMARK 3 B23 (A**2) : .00
REMARK 3
REMARK 3 ESTIMATED OVERALL COORDINATE ERROR.
REMARK 3 ESU BASED ON R VALUE (A) : NULL
REMARK 3 ESU BASED ON FREE R VALUE (A) : .427
REMARK 3 ESU BASED ON MAXIMUM LIKELIHOOD (A) : .327
REMARK 3 ESU FOR B VALUES BASED ON MAXIMUM LIKELIHOOD (A**2) : 39.836

```

REMARK 3
 REMARK 3 CORRELATION COEFFICIENTS.
 REMARK 3 CORRELATION COEFFICIENT FO-FC : .916
 REMARK 3 CORRELATION COEFFICIENT FO-FC FREE : .868
 REMARK 3
 REMARK 3 RMS DEVIATIONS FROM IDEAL VALUES COUNT RMS WEIGHT
 REMARK 3 BOND LENGTHS REFINED ATOMS (A): 8495 ; .011 ; .022
 REMARK 3 BOND LENGTHS OTHERS (A): 5804 ; .001 ; .020
 REMARK 3 BOND ANGLES REFINED ATOMS (DEGREES): 11476 ; 1.279 ; 1.953
 REMARK 3 BOND ANGLES OTHERS (DEGREES): 14093 ; .882 ; 3.000
 REMARK 3 TORSION ANGLES, PERIOD 1 (DEGREES): 1020 ; 7.002 ; 5.000
 REMARK 3 TORSION ANGLES, PERIOD 2 (DEGREES): 435 ; 35.412 ; 24.299
 REMARK 3 TORSION ANGLES, PERIOD 3 (DEGREES): 1525 ; 18.250 ; 15.000
 REMARK 3 TORSION ANGLES, PERIOD 4 (DEGREES): 58 ; 16.811 ; 15.000
 REMARK 3 CHIRAL-CENTER RESTRAINTS (A**3): 1266 ; .070 ; .200
 REMARK 3 GENERAL PLANES REFINED ATOMS (A): 9416 ; .005 ; .020
 REMARK 3 GENERAL PLANES OTHERS (A): 1780 ; .001 ; .020
 REMARK 3
 REMARK 3 ISOTROPIC THERMAL FACTOR RESTRAINTS. COUNT RMS WEIGHT
 REMARK 3 MAIN-CHAIN BOND REFINED ATOMS (A**2): 5104 ; .514 ; 1.500
 REMARK 3 MAIN-CHAIN BOND OTHER ATOMS (A**2): 2068 ; .059 ; 1.500
 REMARK 3 MAIN-CHAIN ANGLE REFINED ATOMS (A**2): 8204 ; 1.000 ; 2.000
 REMARK 3 SIDE-CHAIN BOND REFINED ATOMS (A**2): 3391 ; 1.218 ; 3.000
 REMARK 3 SIDE-CHAIN ANGLE REFINED ATOMS (A**2): 3272 ; 2.157 ; 4.500
 REMARK 3
 REMARK 3 NCS RESTRAINTS STATISTICS
 REMARK 3 NUMBER OF DIFFERENT NCS GROUPS : 1
 REMARK 3
 REMARK 3 NCS GROUP NUMBER : 1
 REMARK 3 CHAIN NAMES : A B
 REMARK 3 NUMBER OF COMPONENTS NCS GROUP : 1
 REMARK 3 COMPONENT C SSSEQI TO C SSSEQI CODE
 REMARK 3 1 A 17 A 541 6
 REMARK 3 1 B 17 B 541 6
 REMARK 3 GROUP CHAIN COUNT RMS WEIGHT
 REMARK 3 LOOSE POSITIONAL 1 1 (A): 7038 ; .37 ; 5.00
 REMARK 3 LOOSE THERMAL 1 1 (A**2): 7038 ; 1.09 ; 10.00
 REMARK 3
 REMARK 3 TWIN DETAILS
 REMARK 3 NUMBER OF TWIN DOMAINS : NULL
 REMARK 3
 REMARK 3
 REMARK 3 TLS DETAILS
 REMARK 3 NUMBER OF TLS GROUPS : 8
 REMARK 3 ATOM RECORD CONTAINS RESIDUAL B FACTORS ONLY
 REMARK 3
 REMARK 3 TLS GROUP : 1
 REMARK 3 NUMBER OF COMPONENTS GROUP : 1
 REMARK 3 COMPONENTS C SSSEQI TO C SSSEQI
 REMARK 3 RESIDUE RANGE : A 17 A 219
 REMARK 3 ORIGIN FOR THE GROUP (A): -64.7667 37.6643 -.0896
 REMARK 3 T TENSOR
 REMARK 3 T11: .0648 T22: .0357
 REMARK 3 T33: .0787 T12: .0200
 REMARK 3 T13: .0129 T23: -.0089
 REMARK 3 L TENSOR
 REMARK 3 L11: 3.7204 L22: 1.5111
 REMARK 3 L33: 2.6701 L12: .5715
 REMARK 3 L13: .6692 L23: -.9699
 REMARK 3 S TENSOR

REMARK 3 S11: .0562 S12: .0478 S13: -.1976
 REMARK 3 S21: -.1702 S22: -.0055 S23: .1376
 REMARK 3 S31: .0900 S32: -.2188 S33: -.0507
 REMARK 3
 REMARK 3 TLS GROUP : 2
 REMARK 3 NUMBER OF COMPONENTS GROUP : 1
 REMARK 3 COMPONENTS C SSSEI TO C SSSEI
 REMARK 3 RESIDUE RANGE : A 220 A 287
 REMARK 3 ORIGIN FOR THE GROUP (A): -59.5787 8.4529 -.7693
 REMARK 3 T TENSOR
 REMARK 3 T11: .1615 T22: .0645
 REMARK 3 T33: .1539 T12: -.0314
 REMARK 3 T13: -.0461 T23: .0198
 REMARK 3 L TENSOR
 REMARK 3 L11: 2.4192 L22: 4.6709
 REMARK 3 L33: .7709 L12: -3.2943
 REMARK 3 L13: -.1814 L23: -.0705
 REMARK 3 S TENSOR
 REMARK 3 S11: .0055 S12: -.0699 S13: -.2073
 REMARK 3 S21: -.1805 S22: .0996 S23: .3781
 REMARK 3 S31: .2596 S32: .0887 S33: -.1051
 REMARK 3
 REMARK 3 TLS GROUP : 3
 REMARK 3 NUMBER OF COMPONENTS GROUP : 1
 REMARK 3 COMPONENTS C SSSEI TO C SSSEI
 REMARK 3 RESIDUE RANGE : A 288 A 374
 REMARK 3 ORIGIN FOR THE GROUP (A): -40.1866 1.6932 .5805
 REMARK 3 T TENSOR
 REMARK 3 T11: .1149 T22: .1003
 REMARK 3 T33: .1629 T12: .0153
 REMARK 3 T13: .0224 T23: .0164
 REMARK 3 L TENSOR
 REMARK 3 L11: .2271 L22: .7399
 REMARK 3 L33: 4.8529 L12: .3413
 REMARK 3 L13: .4755 L23: -.1746
 REMARK 3 S TENSOR
 REMARK 3 S11: -.0449 S12: -.0288 S13: -.1131
 REMARK 3 S21: -.1346 S22: -.0665 S23: -.2749
 REMARK 3 S31: -.0040 S32: .1558 S33: .1114
 REMARK 3
 REMARK 3 TLS GROUP : 4
 REMARK 3 NUMBER OF COMPONENTS GROUP : 1
 REMARK 3 COMPONENTS C SSSEI TO C SSSEI
 REMARK 3 RESIDUE RANGE : A 375 A 541
 REMARK 3 ORIGIN FOR THE GROUP (A): -47.2220 21.5399 6.9217
 REMARK 3 T TENSOR
 REMARK 3 T11: .1551 T22: .1194
 REMARK 3 T33: .1485 T12: -.0666
 REMARK 3 T13: .0275 T23: .0272
 REMARK 3 L TENSOR
 REMARK 3 L11: 2.2352 L22: 2.1698
 REMARK 3 L33: 2.3370 L12: -.4501
 REMARK 3 L13: 2.2662 L23: -.1852
 REMARK 3 S TENSOR
 REMARK 3 S11: .0233 S12: -.3041 S13: -.0323
 REMARK 3 S21: .3375 S22: .0236 S23: .0121
 REMARK 3 S31: .0592 S32: -.2979 S33: -.0469
 REMARK 3
 REMARK 3 TLS GROUP : 5
 REMARK 3 NUMBER OF COMPONENTS GROUP : 1

REMARK 3 COMPONENTS C SSSEI TO C SSSEI
 REMARK 3 RESIDUE RANGE : B 17 B 219
 REMARK 3 ORIGIN FOR THE GROUP (A): -73.9834 -39.9016 -18.5783
 REMARK 3 T TENSOR
 REMARK 3 T11: .0658 T22: .1153
 REMARK 3 T33: .1251 T12: -.0621
 REMARK 3 T13: -.0164 T23: -.0098
 REMARK 3 L TENSOR
 REMARK 3 L11: 4.6230 L22: 1.7260
 REMARK 3 L33: 3.8816 L12: -.4202
 REMARK 3 L13: -1.8646 L23: -.9046
 REMARK 3 S TENSOR
 REMARK 3 S11: -.0685 S12: .0375 S13: -.0003
 REMARK 3 S21: .1931 S22: .0510 S23: -.0097
 REMARK 3 S31: .0317 S32: -.2047 S33: .0175
 REMARK 3
 REMARK 3 TLS GROUP : 6
 REMARK 3 NUMBER OF COMPONENTS GROUP : 1
 REMARK 3 COMPONENTS C SSSEI TO C SSSEI
 REMARK 3 RESIDUE RANGE : B 220 B 287
 REMARK 3 ORIGIN FOR THE GROUP (A): -62.1586 -12.7634 -18.1912
 REMARK 3 T TENSOR
 REMARK 3 T11: .1825 T22: .0804
 REMARK 3 T33: .1512 T12: .0549
 REMARK 3 T13: .0773 T23: .0208
 REMARK 3 L TENSOR
 REMARK 3 L11: 5.4421 L22: 4.0606
 REMARK 3 L33: 1.5369 L12: 4.6706
 REMARK 3 L13: -2.0058 L23: -1.5537
 REMARK 3 S TENSOR
 REMARK 3 S11: .1622 S12: .0431 S13: .3257
 REMARK 3 S21: .1755 S22: .0292 S23: .2977
 REMARK 3 S31: -.1910 S32: -.0506 S33: -.1914
 REMARK 3
 REMARK 3 TLS GROUP : 7
 REMARK 3 NUMBER OF COMPONENTS GROUP : 1
 REMARK 3 COMPONENTS C SSSEI TO C SSSEI
 REMARK 3 RESIDUE RANGE : B 288 B 374
 REMARK 3 ORIGIN FOR THE GROUP (A): -41.6930 -10.8250 -19.6636
 REMARK 3 T TENSOR
 REMARK 3 T11: .1424 T22: .0604
 REMARK 3 T33: .1153 T12: .0184
 REMARK 3 T13: .0146 T23: .0276
 REMARK 3 L TENSOR
 REMARK 3 L11: .6426 L22: .8163
 REMARK 3 L33: 2.3437 L12: -.1831
 REMARK 3 L13: -.5246 L23: .4917
 REMARK 3 S TENSOR
 REMARK 3 S11: .0592 S12: -.0206 S13: .0071
 REMARK 3 S21: .0906 S22: -.0229 S23: -.1585
 REMARK 3 S31: -.0355 S32: .0262 S33: -.0363
 REMARK 3
 REMARK 3 TLS GROUP : 8
 REMARK 3 NUMBER OF COMPONENTS GROUP : 1
 REMARK 3 COMPONENTS C SSSEI TO C SSSEI
 REMARK 3 RESIDUE RANGE : B 375 B 541
 REMARK 3 ORIGIN FOR THE GROUP (A): -53.4886 -28.3212 -26.0670
 REMARK 3 T TENSOR
 REMARK 3 T11: .1107 T22: .1220
 REMARK 3 T33: .1514 T12: .0692

REMARK 3 T13: -.0073 T23: .0518
 REMARK 3 L TENSOR
 REMARK 3 L11: 2.6766 L22: 1.8433
 REMARK 3 L33: 2.6389 L12: .1130
 REMARK 3 L13: -2.4696 L23: .6986
 REMARK 3 S TENSOR
 REMARK 3 S11: .1115 S12: .3882 S13: .0569
 REMARK 3 S21: -.0725 S22: -.0724 S23: .1450
 REMARK 3 S31: -.1453 S32: -.4044 S33: -.0392
 REMARK 3
 REMARK 3
 REMARK 3 BULK SOLVENT MODELLING.
 REMARK 3 METHOD USED : MASK
 REMARK 3 PARAMETERS FOR MASK CALCULATION
 REMARK 3 VDW PROBE RADIUS : 1.40
 REMARK 3 ION PROBE RADIUS : .80
 REMARK 3 SHRINKAGE RADIUS : .80
 REMARK 3
 REMARK 3 OTHER REFINEMENT REMARKS:
 REMARK 3 HYDROGENS HAVE BEEN ADDED IN THE RIDING POSITIONS
 REMARK 3 U VALUES : RESIDUAL ONLY
 REMARK 3
 CISPEP 1 ALA A 446 ARG A 447 .00
 CISPEP 2 GLY A 515 ASP A 516 .00
 CISPEP 3 THR B 22 ASP B 23 .00
 CISPEP 4 ALA B 446 ARG B 447 .00
 CISPEP 5 GLY B 515 ASP B 516 .00
 CRYST1 155.800 155.800 143.690 90.00 90.00 90.00 P 43 21 2
 SCALE1 .006418 .000000 .000000 .000000
 SCALE2 .000000 .006418 .000000 .000000
 SCALE3 .000000 .000000 .006959 .000000
 ATOM 1 N LEU A 17 -63.930 24.416 -19.202 1.00 30.90 N
 ATOM 2 CA LEU A 17 -64.132 23.019 -19.731 1.00 31.43 C
 ATOM 4 CB LEU A 17 -63.308 22.800 -21.021 1.00 31.35 C
 ATOM 7 CG LEU A 17 -64.002 23.016 -22.386 1.00 31.95 C
 ATOM 9 CD1 LEU A 17 -62.989 23.068 -23.550 1.00 31.88 C
 ATOM 13 CD2 LEU A 17 -65.052 21.915 -22.676 1.00 31.95 C
 ATOM 17 C LEU A 17 -63.783 21.952 -18.660 1.00 31.39 C
 ATOM 18 O LEU A 17 -63.142 22.291 -17.651 1.00 31.61 O
 ATOM 22 N LEU A 18 -64.216 20.693 -18.877 1.00 31.11 N
 ATOM 23 CA LEU A 18 -63.986 19.529 -17.957 1.00 30.78 C
 ATOM 25 CB LEU A 18 -62.503 19.072 -17.935 1.00 30.97 C
 ATOM 28 CG LEU A 18 -61.303 19.949 -17.496 1.00 31.10 C
 ATOM 30 CD1 LEU A 18 -61.406 20.430 -16.064 1.00 30.71 C
 ATOM 34 CD2 LEU A 18 -59.977 19.173 -17.703 1.00 31.87 C
 ATOM 38 C LEU A 18 -64.531 19.665 -16.522 1.00 30.38 C
 ATOM 39 O LEU A 18 -64.837 18.667 -15.883 1.00 29.86 O
 ATOM 41 N SER A 19 -64.587 20.900 -16.023 1.00 30.45 N
 ATOM 42 CA SER A 19 -65.335 21.276 -14.815 1.00 30.37 C
 ATOM 44 CB SER A 19 -64.604 22.417 -14.026 1.00 30.02 C
 ATOM 47 OG SER A 19 -64.881 23.741 -14.489 1.00 28.00 O
 ATOM 49 C SER A 19 -66.784 21.654 -15.218 1.00 31.03 C
 ATOM 50 O SER A 19 -67.666 21.711 -14.367 1.00 30.95 O
 ATOM 52 N SER A 20 -67.023 21.880 -16.519 1.00 31.85 N
 ATOM 53 CA SER A 20 -68.355 22.247 -17.051 1.00 32.37 C
 ATOM 55 CB SER A 20 -68.291 22.507 -18.579 1.00 32.39 C
 ATOM 58 OG SER A 20 -67.387 23.542 -18.931 1.00 31.77 O
 ATOM 60 C SER A 20 -69.357 21.124 -16.744 1.00 33.04 C
 ATOM 61 O SER A 20 -69.076 20.254 -15.922 1.00 33.08 O
 ATOM 63 N ASP A 21 -70.522 21.136 -17.393 1.00 33.88 N

ATOM	64	CA	ASP	A	21	-71.512	20.069	-17.206	1.00	34.48	C
ATOM	66	CB	ASP	A	21	-72.907	20.655	-17.127	1.00	34.59	C
ATOM	69	CG	ASP	A	21	-73.022	21.661	-16.006	1.00	35.54	C
ATOM	70	OD1	ASP	A	21	-72.041	21.836	-15.251	1.00	35.31	O
ATOM	71	OD2	ASP	A	21	-74.082	22.289	-15.870	1.00	38.72	O
ATOM	72	C	ASP	A	21	-71.409	18.975	-18.260	1.00	34.94	C
ATOM	73	O	ASP	A	21	-72.134	18.947	-19.269	1.00	34.58	O
ATOM	75	N	THR	A	22	-70.457	18.086	-17.989	1.00	35.66	N
ATOM	76	CA	THR	A	22	-70.340	16.799	-18.657	1.00	36.27	C
ATOM	78	CB	THR	A	22	-68.895	16.576	-19.168	1.00	36.30	C
ATOM	80	OG1	THR	A	22	-67.968	17.278	-18.322	1.00	36.04	O
ATOM	82	CG2	THR	A	22	-68.755	17.088	-20.604	1.00	36.11	C
ATOM	86	C	THR	A	22	-70.792	15.714	-17.648	1.00	36.83	C
ATOM	87	O	THR	A	22	-69.968	15.022	-17.035	1.00	37.05	O
ATOM	89	N	ASP	A	23	-72.121	15.599	-17.494	1.00	37.27	N
ATOM	90	CA	ASP	A	23	-72.790	14.802	-16.441	1.00	37.25	C
ATOM	92	CB	ASP	A	23	-72.962	15.659	-15.167	1.00	37.17	C
ATOM	95	CG	ASP	A	23	-71.625	16.081	-14.549	1.00	37.74	C
ATOM	96	OD1	ASP	A	23	-70.714	15.241	-14.436	1.00	38.88	O
ATOM	97	OD2	ASP	A	23	-71.472	17.256	-14.164	1.00	38.91	O
ATOM	98	C	ASP	A	23	-74.172	14.351	-16.956	1.00	37.18	C
ATOM	99	O	ASP	A	23	-75.137	15.112	-16.846	1.00	37.10	O
ATOM	101	N	GLU	A	24	-74.280	13.134	-17.501	1.00	37.23	N
ATOM	102	CA	GLU	A	24	-75.406	12.815	-18.410	1.00	37.58	C
ATOM	104	CB	GLU	A	24	-74.941	13.007	-19.865	1.00	37.77	C
ATOM	107	CG	GLU	A	24	-74.424	14.417	-20.239	1.00	37.98	C
ATOM	110	CD	GLU	A	24	-74.121	14.556	-21.746	1.00	38.35	C
ATOM	111	OE1	GLU	A	24	-74.284	13.564	-22.503	1.00	38.12	O
ATOM	112	OE2	GLU	A	24	-73.721	15.661	-22.174	1.00	38.18	O
ATOM	113	C	GLU	A	24	-76.139	11.440	-18.323	1.00	37.79	C
ATOM	114	O	GLU	A	24	-77.323	11.397	-17.960	1.00	37.90	O
ATOM	116	N	SER	A	25	-75.462	10.345	-18.696	1.00	37.92	N
ATOM	117	CA	SER	A	25	-76.140	9.070	-19.048	1.00	38.09	C
ATOM	119	CB	SER	A	25	-75.167	8.076	-19.723	1.00	38.11	C
ATOM	122	OG	SER	A	25	-74.503	7.239	-18.787	1.00	38.02	O
ATOM	124	C	SER	A	25	-76.917	8.374	-17.907	1.00	38.42	C
ATOM	125	O	SER	A	25	-76.635	8.581	-16.724	1.00	38.15	O
ATOM	127	N	ILE	A	26	-77.861	7.511	-18.312	1.00	39.00	N
ATOM	128	CA	ILE	A	26	-78.993	7.042	-17.476	1.00	39.30	C
ATOM	130	CB	ILE	A	26	-78.594	5.965	-16.429	1.00	39.34	C
ATOM	132	CG1	ILE	A	26	-77.892	4.791	-17.120	1.00	39.37	C
ATOM	135	CD1	ILE	A	26	-77.766	3.531	-16.253	1.00	39.46	C
ATOM	139	CG2	ILE	A	26	-79.836	5.423	-15.714	1.00	39.29	C
ATOM	143	C	ILE	A	26	-79.716	8.260	-16.852	1.00	39.66	C
ATOM	144	O	ILE	A	26	-79.274	8.831	-15.838	1.00	39.49	O
ATOM	146	N	GLU	A	27	-80.837	8.623	-17.486	1.00	39.99	N
ATOM	147	CA	GLU	A	27	-81.450	9.965	-17.383	1.00	40.11	C
ATOM	149	CB	GLU	A	27	-82.395	10.183	-18.587	1.00	40.19	C
ATOM	152	CG	GLU	A	27	-81.632	10.311	-19.904	1.00	40.60	C
ATOM	155	CD	GLU	A	27	-82.528	10.434	-21.116	1.00	41.06	C
ATOM	156	OE1	GLU	A	27	-83.328	9.502	-21.367	1.00	41.10	O
ATOM	157	OE2	GLU	A	27	-82.409	11.458	-21.830	1.00	41.36	O
ATOM	158	C	GLU	A	27	-82.144	10.264	-16.042	1.00	40.02	C
ATOM	159	O	GLU	A	27	-81.977	9.518	-15.071	1.00	40.01	O
ATOM	161	N	VAL	A	28	-82.899	11.370	-16.005	1.00	39.84	N
ATOM	162	CA	VAL	A	28	-83.463	11.960	-14.779	1.00	39.65	C
ATOM	164	CB	VAL	A	28	-83.884	10.905	-13.694	1.00	39.78	C
ATOM	166	CG1	VAL	A	28	-84.472	11.596	-12.456	1.00	39.87	C
ATOM	170	CG2	VAL	A	28	-84.887	9.880	-14.274	1.00	39.49	C
ATOM	174	C	VAL	A	28	-82.469	12.980	-14.206	1.00	39.46	C

ATOM	175	O	VAL A	28	-82.798	13.717	-13.270	1.00	39.38	O
ATOM	177	N	HIS A	29	-81.264	13.025	-14.786	1.00	39.30	N
ATOM	178	CA	HIS A	29	-80.236	14.006	-14.412	1.00	39.23	C
ATOM	180	CB	HIS A	29	-78.866	13.344	-14.186	1.00	39.45	C
ATOM	183	CG	HIS A	29	-78.910	12.139	-13.296	1.00	41.29	C
ATOM	184	ND1	HIS A	29	-78.407	12.140	-12.007	1.00	42.51	N
ATOM	186	CE1	HIS A	29	-78.585	10.942	-11.472	1.00	43.37	C
ATOM	188	NE2	HIS A	29	-79.181	10.164	-12.365	1.00	43.18	N
ATOM	190	CD2	HIS A	29	-79.393	10.888	-13.515	1.00	42.55	C
ATOM	192	C	HIS A	29	-80.090	15.110	-15.464	1.00	38.60	C
ATOM	193	O	HIS A	29	-79.049	15.774	-15.499	1.00	38.61	O
ATOM	195	N	LYS A	30	-81.110	15.320	-16.311	1.00	37.87	N
ATOM	196	CA	LYS A	30	-81.265	16.610	-17.001	1.00	37.22	C
ATOM	198	CB	LYS A	30	-82.283	16.542	-18.140	1.00	37.09	C
ATOM	201	CG	LYS A	30	-81.764	15.794	-19.360	1.00	37.00	C
ATOM	204	CD	LYS A	30	-82.831	15.653	-20.457	1.00	36.89	C
ATOM	207	CE	LYS A	30	-82.413	14.637	-21.531	1.00	36.48	C
ATOM	210	NZ	LYS A	30	-83.422	14.458	-22.611	1.00	35.46	N
ATOM	214	C	LYS A	30	-81.648	17.649	-15.938	1.00	36.88	C
ATOM	215	O	LYS A	30	-82.516	18.499	-16.136	1.00	36.71	O
ATOM	217	N	ASP A	31	-80.982	17.506	-14.788	1.00	36.58	N
ATOM	218	CA	ASP A	31	-80.895	18.476	-13.720	1.00	36.18	C
ATOM	220	CB	ASP A	31	-80.518	17.751	-12.396	1.00	36.20	C
ATOM	223	CG	ASP A	31	-79.863	18.669	-11.342	1.00	36.65	C
ATOM	224	OD1	ASP A	31	-80.055	19.901	-11.390	1.00	38.17	O
ATOM	225	OD2	ASP A	31	-79.147	18.150	-10.449	1.00	35.46	O
ATOM	226	C	ASP A	31	-79.811	19.426	-14.226	1.00	35.74	C
ATOM	227	O	ASP A	31	-78.628	19.279	-13.918	1.00	35.44	O
ATOM	229	N	LYS A	32	-80.220	20.343	-15.096	1.00	35.30	N
ATOM	230	CA	LYS A	32	-79.360	21.449	-15.507	1.00	34.69	C
ATOM	232	CB	LYS A	32	-78.698	21.188	-16.849	1.00	34.63	C
ATOM	235	CG	LYS A	32	-77.699	20.042	-16.765	1.00	34.51	C
ATOM	238	CD	LYS A	32	-76.953	19.853	-18.078	1.00	34.48	C
ATOM	241	CE	LYS A	32	-76.859	18.387	-18.503	1.00	33.77	C
ATOM	244	NZ	LYS A	32	-77.002	18.244	-19.985	1.00	32.92	N
ATOM	248	C	LYS A	32	-80.162	22.745	-15.459	1.00	34.01	C
ATOM	249	O	LYS A	32	-80.468	23.380	-16.473	1.00	33.41	O
ATOM	251	N	ALA A	33	-80.540	23.051	-14.219	1.00	33.21	N
ATOM	252	CA	ALA A	33	-80.785	24.381	-13.766	1.00	32.60	C
ATOM	254	CB	ALA A	33	-81.619	24.347	-12.484	1.00	32.28	C
ATOM	258	C	ALA A	33	-79.411	25.031	-13.524	1.00	32.31	C
ATOM	259	O	ALA A	33	-79.335	26.081	-12.901	1.00	32.80	O
ATOM	261	N	LYS A	34	-78.323	24.389	-13.968	1.00	31.70	N
ATOM	262	CA	LYS A	34	-77.035	25.071	-14.186	1.00	31.12	C
ATOM	264	CB	LYS A	34	-75.863	24.089	-14.128	1.00	31.23	C
ATOM	267	CG	LYS A	34	-75.118	24.077	-12.791	1.00	31.76	C
ATOM	270	CD	LYS A	34	-73.584	24.412	-12.921	1.00	31.49	C
ATOM	273	CE	LYS A	34	-72.685	23.269	-12.506	1.00	30.50	C
ATOM	276	NZ	LYS A	34	-72.951	22.090	-13.348	1.00	29.53	N
ATOM	280	C	LYS A	34	-77.024	25.884	-15.518	1.00	30.50	C
ATOM	281	O	LYS A	34	-76.977	25.353	-16.625	1.00	29.82	O
ATOM	283	N	LYS A	35	-76.927	27.191	-15.337	1.00	29.97	N
ATOM	284	CA	LYS A	35	-77.632	28.233	-16.098	1.00	29.41	C
ATOM	286	CB	LYS A	35	-78.852	27.704	-16.867	1.00	29.53	C
ATOM	289	CG	LYS A	35	-80.166	27.594	-16.077	1.00	30.18	C
ATOM	292	CD	LYS A	35	-81.030	28.858	-16.203	1.00	31.25	C
ATOM	295	CE	LYS A	35	-82.375	28.721	-15.486	1.00	31.85	C
ATOM	298	NZ	LYS A	35	-83.298	29.865	-15.783	1.00	31.73	N
ATOM	302	C	LYS A	35	-78.053	29.277	-15.026	1.00	28.53	C
ATOM	303	O	LYS A	35	-78.246	30.460	-15.307	1.00	28.53	O

ATOM	305	N	LEU	A	36	-78.225	28.791	-13.796	1.00	27.29	N
ATOM	306	CA	LEU	A	36	-77.998	29.576	-12.594	1.00	26.21	C
ATOM	308	CB	LEU	A	36	-78.088	28.679	-11.364	1.00	25.76	C
ATOM	311	CG	LEU	A	36	-79.468	28.350	-10.832	1.00	24.08	C
ATOM	313	CD1	LEU	A	36	-79.376	27.233	-9.836	1.00	22.16	C
ATOM	317	CD2	LEU	A	36	-80.051	29.586	-10.213	1.00	23.45	C
ATOM	321	C	LEU	A	36	-76.587	30.137	-12.665	1.00	25.93	C
ATOM	322	O	LEU	A	36	-76.290	31.221	-12.137	1.00	26.11	O
ATOM	324	N	GLU	A	37	-75.714	29.332	-13.260	1.00	25.24	N
ATOM	325	CA	GLU	A	37	-74.381	29.739	-13.650	1.00	24.85	C
ATOM	327	CB	GLU	A	37	-73.678	28.566	-14.313	1.00	24.74	C
ATOM	330	CG	GLU	A	37	-72.182	28.667	-14.285	1.00	23.93	C
ATOM	333	CD	GLU	A	37	-71.519	27.490	-14.937	1.00	22.32	C
ATOM	334	OE1	GLU	A	37	-72.203	26.738	-15.657	1.00	20.57	O
ATOM	335	OE2	GLU	A	37	-70.303	27.322	-14.728	1.00	22.27	O
ATOM	336	C	GLU	A	37	-74.391	30.927	-14.600	1.00	24.69	C
ATOM	337	O	GLU	A	37	-73.666	31.891	-14.389	1.00	24.49	O
ATOM	339	N	ALA	A	38	-75.203	30.850	-15.650	1.00	24.68	N
ATOM	340	CA	ALA	A	38	-75.352	31.960	-16.599	1.00	24.87	C
ATOM	342	CB	ALA	A	38	-76.453	31.656	-17.609	1.00	24.53	C
ATOM	346	C	ALA	A	38	-75.643	33.278	-15.880	1.00	25.08	C
ATOM	347	O	ALA	A	38	-75.009	34.295	-16.135	1.00	25.21	O
ATOM	349	N	GLU	A	39	-76.591	33.228	-14.956	1.00	25.35	N
ATOM	350	CA	GLU	A	39	-77.050	34.400	-14.219	1.00	25.51	C
ATOM	352	CB	GLU	A	39	-78.283	34.014	-13.396	1.00	25.95	C
ATOM	355	CG	GLU	A	39	-79.302	35.119	-13.136	1.00	27.03	C
ATOM	358	CD	GLU	A	39	-80.715	34.557	-12.906	1.00	28.70	C
ATOM	359	OE1	GLU	A	39	-81.006	33.416	-13.346	1.00	27.40	O
ATOM	360	OE2	GLU	A	39	-81.540	35.266	-12.289	1.00	31.12	O
ATOM	361	C	GLU	A	39	-75.970	34.966	-13.306	1.00	25.11	C
ATOM	362	O	GLU	A	39	-75.870	36.167	-13.164	1.00	25.28	O
ATOM	364	N	VAL	A	40	-75.182	34.107	-12.672	1.00	24.86	N
ATOM	365	CA	VAL	A	40	-74.079	34.568	-11.824	1.00	24.69	C
ATOM	367	CB	VAL	A	40	-73.511	33.420	-10.952	1.00	24.56	C
ATOM	369	CG1	VAL	A	40	-72.239	33.857	-10.240	1.00	23.69	C
ATOM	373	CG2	VAL	A	40	-74.553	32.955	-9.955	1.00	24.68	C
ATOM	377	C	VAL	A	40	-72.971	35.148	-12.698	1.00	24.76	C
ATOM	378	O	VAL	A	40	-72.337	36.139	-12.362	1.00	24.16	O
ATOM	380	N	ARG	A	41	-72.744	34.506	-13.831	1.00	25.25	N
ATOM	381	CA	ARG	A	41	-71.727	34.948	-14.769	1.00	25.66	C
ATOM	383	CB	ARG	A	41	-71.576	33.927	-15.896	1.00	25.88	C
ATOM	386	CG	ARG	A	41	-70.726	34.385	-17.062	1.00	27.04	C
ATOM	389	CD	ARG	A	41	-71.519	35.178	-18.095	1.00	27.90	C
ATOM	392	NE	ARG	A	41	-70.653	35.635	-19.180	1.00	29.20	N
ATOM	394	CZ	ARG	A	41	-70.946	36.615	-20.034	1.00	29.94	C
ATOM	395	NH1	ARG	A	41	-70.077	36.945	-20.980	1.00	29.86	N
ATOM	398	NH2	ARG	A	41	-72.096	37.273	-19.957	1.00	30.54	N
ATOM	401	C	ARG	A	41	-72.104	36.301	-15.335	1.00	25.60	C
ATOM	402	O	ARG	A	41	-71.237	37.113	-15.612	1.00	25.77	O
ATOM	404	N	ARG	A	42	-73.400	36.521	-15.537	1.00	25.50	N
ATOM	405	CA	ARG	A	42	-73.900	37.810	-15.979	1.00	25.39	C
ATOM	407	CB	ARG	A	42	-75.389	37.723	-16.313	1.00	25.09	C
ATOM	410	CG	ARG	A	42	-76.077	39.062	-16.445	1.00	23.84	C
ATOM	413	CD	ARG	A	42	-77.459	38.906	-16.980	1.00	22.37	C
ATOM	416	NE	ARG	A	42	-78.448	38.554	-15.965	1.00	21.69	N
ATOM	418	CZ	ARG	A	42	-79.705	38.206	-16.248	1.00	22.99	C
ATOM	419	NH1	ARG	A	42	-80.122	38.146	-17.517	1.00	24.16	N
ATOM	422	NH2	ARG	A	42	-80.557	37.900	-15.275	1.00	22.94	N
ATOM	425	C	ARG	A	42	-73.666	38.889	-14.925	1.00	26.00	C
ATOM	426	O	ARG	A	42	-73.173	39.959	-15.244	1.00	26.17	O

ATOM	428	N	GLU A 43	-74.014	38.621	-13.673	1.00	26.66	N
ATOM	429	CA	GLU A 43	-73.930	39.662	-12.650	1.00	27.47	C
ATOM	431	CB	GLU A 43	-74.749	39.304	-11.401	1.00	27.97	C
ATOM	434	CG	GLU A 43	-76.274	39.266	-11.661	1.00	31.05	C
ATOM	437	CD	GLU A 43	-76.869	40.636	-12.063	1.00	34.68	C
ATOM	438	OE1	GLU A 43	-76.804	41.571	-11.222	1.00	37.63	O
ATOM	439	OE2	GLU A 43	-77.398	40.773	-13.205	1.00	34.52	O
ATOM	440	C	GLU A 43	-72.494	40.035	-12.275	1.00	27.15	C
ATOM	441	O	GLU A 43	-72.292	41.118	-11.716	1.00	27.49	O
ATOM	443	N	ILE A 44	-71.517	39.166	-12.588	1.00	26.57	N
ATOM	444	CA	ILE A 44	-70.090	39.456	-12.334	1.00	25.80	C
ATOM	446	CB	ILE A 44	-69.206	38.187	-12.187	1.00	25.52	C
ATOM	448	CG1	ILE A 44	-69.624	37.320	-11.010	1.00	24.69	C
ATOM	451	CD1	ILE A 44	-68.828	36.061	-10.899	1.00	23.08	C
ATOM	455	CG2	ILE A 44	-67.790	38.581	-11.916	1.00	25.73	C
ATOM	459	C	ILE A 44	-69.522	40.286	-13.472	1.00	25.40	C
ATOM	460	O	ILE A 44	-68.745	41.211	-13.236	1.00	25.26	O
ATOM	462	N	ASN A 45	-69.912	39.945	-14.700	1.00	25.14	N
ATOM	463	CA	ASN A 45	-69.480	40.669	-15.906	1.00	25.02	C
ATOM	465	CB	ASN A 45	-69.696	39.813	-17.151	1.00	24.69	C
ATOM	468	CG	ASN A 45	-68.662	38.735	-17.302	1.00	24.01	C
ATOM	469	OD1	ASN A 45	-67.470	39.013	-17.395	1.00	23.81	O
ATOM	470	ND2	ASN A 45	-69.111	37.490	-17.344	1.00	23.28	N
ATOM	473	C	ASN A 45	-70.177	42.014	-16.114	1.00	25.52	C
ATOM	474	O	ASN A 45	-69.704	42.836	-16.890	1.00	25.40	O
ATOM	476	N	ASN A 46	-71.308	42.210	-15.437	1.00	26.33	N
ATOM	477	CA	ASN A 46	-72.085	43.453	-15.468	1.00	26.72	C
ATOM	479	CB	ASN A 46	-73.102	43.404	-14.322	1.00	26.31	C
ATOM	482	CG	ASN A 46	-73.935	44.641	-14.209	1.00	24.24	C
ATOM	483	OD1	ASN A 46	-73.948	45.472	-15.091	1.00	22.38	O
ATOM	484	ND2	ASN A 46	-74.645	44.766	-13.106	1.00	21.64	N
ATOM	487	C	ASN A 46	-71.185	44.680	-15.349	1.00	28.11	C
ATOM	488	O	ASN A 46	-70.603	44.933	-14.304	1.00	28.16	O
ATOM	490	N	GLU A 47	-71.084	45.441	-16.433	1.00	29.93	N
ATOM	491	CA	GLU A 47	-70.128	46.558	-16.544	1.00	31.30	C
ATOM	493	CB	GLU A 47	-69.933	46.969	-18.013	1.00	31.28	C
ATOM	496	CG	GLU A 47	-69.737	45.795	-19.020	1.00	32.37	C
ATOM	499	CD	GLU A 47	-71.060	45.127	-19.545	1.00	32.69	C
ATOM	500	OE1	GLU A 47	-72.055	45.833	-19.816	1.00	32.16	O
ATOM	501	OE2	GLU A 47	-71.087	43.880	-19.699	1.00	32.41	O
ATOM	502	C	GLU A 47	-70.582	47.774	-15.729	1.00	32.47	C
ATOM	503	O	GLU A 47	-69.760	48.525	-15.230	1.00	32.86	O
ATOM	505	N	LYS A 48	-71.898	47.945	-15.585	1.00	33.93	N
ATOM	506	CA	LYS A 48	-72.491	49.075	-14.861	1.00	34.85	C
ATOM	508	CB	LYS A 48	-73.660	49.672	-15.684	1.00	35.08	C
ATOM	511	CG	LYS A 48	-73.221	50.119	-17.126	1.00	36.58	C
ATOM	514	CD	LYS A 48	-74.035	51.277	-17.751	1.00	37.93	C
ATOM	517	CE	LYS A 48	-75.353	50.799	-18.415	1.00	39.02	C
ATOM	520	NZ	LYS A 48	-75.159	49.990	-19.669	1.00	39.11	N
ATOM	524	C	LYS A 48	-72.926	48.632	-13.465	1.00	35.26	C
ATOM	525	O	LYS A 48	-74.011	48.954	-13.015	1.00	35.04	O
ATOM	527	N	ALA A 49	-72.051	47.890	-12.790	1.00	36.29	N
ATOM	528	CA	ALA A 49	-72.321	47.348	-11.456	1.00	37.18	C
ATOM	530	CB	ALA A 49	-71.828	45.923	-11.353	1.00	37.20	C
ATOM	534	C	ALA A 49	-71.618	48.191	-10.418	1.00	37.88	C
ATOM	535	O	ALA A 49	-70.480	48.610	-10.638	1.00	38.17	O
ATOM	537	N	GLU A 50	-72.278	48.404	-9.280	1.00	38.53	N
ATOM	538	CA	GLU A 50	-71.730	49.242	-8.213	1.00	39.08	C
ATOM	540	CB	GLU A 50	-72.790	49.545	-7.146	1.00	39.43	C
ATOM	543	CG	GLU A 50	-72.708	50.953	-6.600	1.00	40.86	C

ATOM	546	CD	GLU	A	50	-73.089	52.010	-7.639	1.00	42.77	C
ATOM	547	OE1	GLU	A	50	-74.225	52.530	-7.570	1.00	45.07	O
ATOM	548	OE2	GLU	A	50	-72.265	52.320	-8.530	1.00	43.35	O
ATOM	549	C	GLU	A	50	-70.560	48.521	-7.600	1.00	38.89	C
ATOM	550	O	GLU	A	50	-70.699	47.378	-7.201	1.00	39.05	O
ATOM	552	N	PHE	A	51	-69.411	49.181	-7.536	1.00	38.99	N
ATOM	553	CA	PHE	A	51	-68.161	48.502	-7.196	1.00	39.52	C
ATOM	555	CB	PHE	A	51	-67.000	49.494	-7.098	1.00	40.05	C
ATOM	558	CG	PHE	A	51	-66.460	49.948	-8.431	1.00	42.80	C
ATOM	559	CD1	PHE	A	51	-66.135	49.014	-9.435	1.00	45.21	C
ATOM	561	CE1	PHE	A	51	-65.612	49.425	-10.675	1.00	46.12	C
ATOM	563	CZ	PHE	A	51	-65.403	50.790	-10.917	1.00	47.21	C
ATOM	565	CE2	PHE	A	51	-65.726	51.741	-9.913	1.00	46.72	C
ATOM	567	CD2	PHE	A	51	-66.247	51.310	-8.679	1.00	45.05	C
ATOM	569	C	PHE	A	51	-68.223	47.687	-5.905	1.00	39.05	C
ATOM	570	O	PHE	A	51	-67.897	46.506	-5.900	1.00	39.28	O
ATOM	572	N	LEU	A	52	-68.648	48.307	-4.813	1.00	38.42	N
ATOM	573	CA	LEU	A	52	-68.709	47.610	-3.528	1.00	37.94	C
ATOM	575	CB	LEU	A	52	-69.084	48.594	-2.413	1.00	38.53	C
ATOM	578	CG	LEU	A	52	-68.057	49.726	-2.181	1.00	40.46	C
ATOM	580	CD1	LEU	A	52	-68.691	51.101	-1.787	1.00	41.88	C
ATOM	584	CD2	LEU	A	52	-67.020	49.274	-1.141	1.00	41.83	C
ATOM	588	C	LEU	A	52	-69.667	46.406	-3.532	1.00	36.72	C
ATOM	589	O	LEU	A	52	-69.517	45.494	-2.724	1.00	36.84	O
ATOM	591	N	THR	A	53	-70.649	46.398	-4.431	1.00	35.37	N
ATOM	592	CA	THR	A	53	-71.545	45.247	-4.569	1.00	34.11	C
ATOM	594	CB	THR	A	53	-72.867	45.594	-5.277	1.00	33.87	C
ATOM	596	OG1	THR	A	53	-73.339	46.863	-4.828	1.00	33.49	O
ATOM	598	CG2	THR	A	53	-73.917	44.535	-4.983	1.00	33.70	C
ATOM	602	C	THR	A	53	-70.853	44.148	-5.359	1.00	33.04	C
ATOM	603	O	THR	A	53	-70.893	42.977	-4.973	1.00	32.90	O
ATOM	605	N	LEU	A	54	-70.228	44.528	-6.467	1.00	31.76	N
ATOM	606	CA	LEU	A	54	-69.469	43.578	-7.286	1.00	31.13	C
ATOM	608	CB	LEU	A	54	-68.721	44.308	-8.392	1.00	30.92	C
ATOM	611	CG	LEU	A	54	-68.028	43.437	-9.424	1.00	30.27	C
ATOM	613	CD1	LEU	A	54	-69.034	42.769	-10.310	1.00	29.56	C
ATOM	617	CD2	LEU	A	54	-67.108	44.306	-10.242	1.00	31.15	C
ATOM	621	C	LEU	A	54	-68.466	42.808	-6.443	1.00	30.60	C
ATOM	622	O	LEU	A	54	-68.378	41.587	-6.543	1.00	30.30	O
ATOM	624	N	LEU	A	55	-67.725	43.551	-5.618	1.00	30.09	N
ATOM	625	CA	LEU	A	55	-66.730	43.000	-4.695	1.00	29.47	C
ATOM	627	CB	LEU	A	55	-66.006	44.126	-3.944	1.00	29.44	C
ATOM	630	CG	LEU	A	55	-65.069	45.017	-4.781	1.00	29.74	C
ATOM	632	CD1	LEU	A	55	-64.609	46.255	-4.003	1.00	29.43	C
ATOM	636	CD2	LEU	A	55	-63.859	44.234	-5.286	1.00	29.64	C
ATOM	640	C	LEU	A	55	-67.340	42.024	-3.696	1.00	28.87	C
ATOM	641	O	LEU	A	55	-66.746	40.991	-3.383	1.00	28.98	O
ATOM	643	N	GLU	A	56	-68.525	42.330	-3.198	1.00	28.01	N
ATOM	644	CA	GLU	A	56	-69.160	41.419	-2.256	1.00	27.64	C
ATOM	646	CB	GLU	A	56	-70.172	42.173	-1.395	1.00	28.19	C
ATOM	649	CG	GLU	A	56	-69.497	43.167	-4.459	1.00	29.75	C
ATOM	652	CD	GLU	A	56	-70.457	43.858	.478	1.00	33.04	C
ATOM	653	OE1	GLU	A	56	-71.667	43.538	.481	1.00	34.46	O
ATOM	654	OE2	GLU	A	56	-69.988	44.737	1.226	1.00	36.78	O
ATOM	655	C	GLU	A	56	-69.776	40.200	-2.945	1.00	26.37	C
ATOM	656	O	GLU	A	56	-69.914	39.146	-2.333	1.00	26.02	O
ATOM	658	N	LEU	A	57	-70.134	40.347	-4.218	1.00	25.30	N
ATOM	659	CA	LEU	A	57	-70.569	39.212	-5.036	1.00	24.17	C
ATOM	661	CB	LEU	A	57	-71.125	39.680	-6.382	1.00	23.71	C
ATOM	664	CG	LEU	A	57	-71.417	38.568	-7.390	1.00	22.27	C

ATOM	666	CD1 LEU A 57	-72.515	37.675	-6.858	1.00	19.31	C
ATOM	670	CD2 LEU A 57	-71.752	39.167	-8.768	1.00	20.30	C
ATOM	674	C LEU A 57	-69.404	38.267	-5.284	1.00	23.69	C
ATOM	675	O LEU A 57	-69.547	37.051	-5.136	1.00	23.83	O
ATOM	677	N ILE A 58	-68.261	38.822	-5.682	1.00	22.82	N
ATOM	678	CA ILE A 58	-67.070	38.019	-5.925	1.00	22.39	C
ATOM	680	CB ILE A 58	-65.884	38.899	-6.355	1.00	22.21	C
ATOM	682	CG1 ILE A 58	-66.088	39.386	-7.792	1.00	21.83	C
ATOM	685	CD1 ILE A 58	-65.002	40.314	-8.320	1.00	20.43	C
ATOM	689	CG2 ILE A 58	-64.565	38.131	-6.245	1.00	22.64	C
ATOM	693	C ILE A 58	-66.719	37.227	-4.662	1.00	22.35	C
ATOM	694	O ILE A 58	-66.520	36.002	-4.706	1.00	21.95	O
ATOM	696	N ASP A 59	-66.676	37.951	-3.543	1.00	22.32	N
ATOM	697	CA ASP A 59	-66.396	37.396	-2.215	1.00	22.26	C
ATOM	699	CB ASP A 59	-66.570	38.500	-1.171	1.00	22.67	C
ATOM	702	CG ASP A 59	-66.126	38.092	.213	1.00	23.91	C
ATOM	703	OD1 ASP A 59	-65.271	37.182	.338	1.00	23.77	O
ATOM	704	OD2 ASP A 59	-66.641	38.718	1.178	1.00	27.12	O
ATOM	705	C ASP A 59	-67.304	36.216	-1.901	1.00	21.79	C
ATOM	706	O ASP A 59	-66.828	35.102	-1.668	1.00	21.71	O
ATOM	708	N ASN A 60	-68.610	36.466	-1.926	1.00	21.51	N
ATOM	709	CA ASN A 60	-69.619	35.418	-1.759	1.00	21.28	C
ATOM	711	CB ASN A 60	-71.012	36.006	-1.958	1.00	21.25	C
ATOM	714	CG ASN A 60	-71.476	36.827	-.776	1.00	21.82	C
ATOM	715	OD1 ASN A 60	-71.067	36.598	.364	1.00	20.81	O
ATOM	716	ND2 ASN A 60	-72.361	37.783	-1.043	1.00	23.63	N
ATOM	719	C ASN A 60	-69.445	34.240	-2.725	1.00	21.22	C
ATOM	720	O ASN A 60	-69.444	33.053	-2.303	1.00	20.95	O
ATOM	722	N VAL A 61	-69.308	34.566	-4.018	1.00	20.83	N
ATOM	723	CA VAL A 61	-69.147	33.536	-5.043	1.00	20.57	C
ATOM	725	CB VAL A 61	-68.915	34.122	-6.466	1.00	20.53	C
ATOM	727	CG1 VAL A 61	-68.382	33.042	-7.430	1.00	19.44	C
ATOM	731	CG2 VAL A 61	-70.190	34.758	-7.006	1.00	19.75	C
ATOM	735	C VAL A 61	-67.981	32.646	-4.644	1.00	20.62	C
ATOM	736	O VAL A 61	-68.080	31.431	-4.699	1.00	20.56	O
ATOM	738	N GLN A 62	-66.888	33.255	-4.209	1.00	20.76	N
ATOM	739	CA GLN A 62	-65.703	32.481	-3.898	1.00	21.18	C
ATOM	741	CB GLN A 62	-64.461	33.379	-3.824	1.00	21.39	C
ATOM	744	CG GLN A 62	-64.007	33.920	-5.165	1.00	21.64	C
ATOM	747	CD GLN A 62	-62.608	34.512	-5.135	1.00	22.04	C
ATOM	748	OE1 GLN A 62	-61.890	34.483	-6.145	1.00	23.18	O
ATOM	749	NE2 GLN A 62	-62.217	35.062	-3.989	1.00	20.39	N
ATOM	752	C GLN A 62	-65.891	31.681	-2.607	1.00	21.09	C
ATOM	753	O GLN A 62	-65.544	30.485	-2.554	1.00	21.30	O
ATOM	755	N ARG A 63	-66.448	32.325	-1.583	1.00	20.63	N
ATOM	756	CA ARG A 63	-66.590	31.677	-.273	1.00	20.54	C
ATOM	758	CB ARG A 63	-67.025	32.686	.786	1.00	20.57	C
ATOM	761	CG ARG A 63	-66.031	33.823	.930	1.00	22.55	C
ATOM	764	CD ARG A 63	-66.214	34.632	2.179	1.00	24.93	C
ATOM	767	NE ARG A 63	-66.123	33.781	3.355	1.00	27.75	N
ATOM	769	CZ ARG A 63	-66.496	34.151	4.573	1.00	30.24	C
ATOM	770	NH1 ARG A 63	-66.971	35.385	4.782	1.00	31.79	N
ATOM	773	NH2 ARG A 63	-66.399	33.285	5.581	1.00	30.07	N
ATOM	776	C ARG A 63	-67.557	30.499	-.338	1.00	19.81	C
ATOM	777	O ARG A 63	-67.283	29.431	.209	1.00	19.67	O
ATOM	779	N LEU A 64	-68.661	30.684	-1.057	1.00	19.04	N
ATOM	780	CA LEU A 64	-69.636	29.617	-1.246	1.00	18.44	C
ATOM	782	CB LEU A 64	-70.865	30.167	-1.954	1.00	18.59	C
ATOM	785	CG LEU A 64	-71.662	31.170	-1.140	1.00	18.32	C
ATOM	787	CD1 LEU A 64	-72.669	31.852	-2.034	1.00	18.10	C

ATOM	791	CD2 LEU A 64	-72.334	30.444	.004	1.00	18.29	C
ATOM	795	C LEU A 64	-69.089	28.436	-2.037	1.00	17.80	C
ATOM	796	O LEU A 64	-69.802	27.453	-2.252	1.00	17.79	O
ATOM	798	N GLY A 65	-67.850	28.573	-2.515	1.00	17.07	N
ATOM	799	CA GLY A 65	-67.084	27.480	-3.081	1.00	16.27	C
ATOM	802	C GLY A 65	-67.137	27.406	-4.588	1.00	15.77	C
ATOM	803	O GLY A 65	-66.893	26.347	-5.154	1.00	15.78	O
ATOM	805	N LEU A 66	-67.439	28.521	-5.246	1.00	15.24	N
ATOM	806	CA LEU A 66	-67.621	28.537	-6.701	1.00	14.93	C
ATOM	808	CB LEU A 66	-68.977	29.146	-7.021	1.00	14.84	C
ATOM	811	CG LEU A 66	-70.204	28.334	-6.619	1.00	13.97	C
ATOM	813	CD1 LEU A 66	-71.391	29.260	-6.650	1.00	13.74	C
ATOM	817	CD2 LEU A 66	-70.415	27.125	-7.534	1.00	11.22	C
ATOM	821	C LEU A 66	-66.533	29.314	-7.454	1.00	15.09	C
ATOM	822	O LEU A 66	-66.621	29.497	-8.678	1.00	14.81	O
ATOM	824	N GLY A 67	-65.510	29.758	-6.720	1.00	15.21	N
ATOM	825	CA GLY A 67	-64.404	30.526	-7.280	1.00	14.98	C
ATOM	828	C GLY A 67	-63.873	29.952	-8.569	1.00	14.60	C
ATOM	829	O GLY A 67	-63.891	30.619	-9.579	1.00	14.63	O
ATOM	831	N TYR A 68	-63.419	28.706	-8.523	1.00	14.57	N
ATOM	832	CA TYR A 68	-62.831	28.029	-9.688	1.00	14.72	C
ATOM	834	CB TYR A 68	-62.608	26.539	-9.372	1.00	14.73	C
ATOM	837	CG TYR A 68	-63.858	25.689	-9.330	1.00	12.47	C
ATOM	838	CD1 TYR A 68	-64.163	24.826	-10.366	1.00	10.68	C
ATOM	840	CE1 TYR A 68	-65.310	24.043	-10.340	1.00	10.55	C
ATOM	842	CZ TYR A 68	-66.163	24.119	-9.258	1.00	10.41	C
ATOM	843	OH TYR A 68	-67.311	23.349	-9.212	1.00	6.97	O
ATOM	845	CE2 TYR A 68	-65.864	24.979	-8.210	1.00	11.44	C
ATOM	847	CD2 TYR A 68	-64.722	25.750	-8.252	1.00	11.29	C
ATOM	849	C TYR A 68	-63.648	28.141	-10.980	1.00	15.31	C
ATOM	850	O TYR A 68	-63.106	28.230	-12.083	1.00	15.01	O
ATOM	852	N ARG A 69	-64.959	28.136	-10.816	1.00	16.03	N
ATOM	853	CA ARG A 69	-65.884	28.083	-11.922	1.00	16.67	C
ATOM	855	CB ARG A 69	-67.224	27.649	-11.348	1.00	16.58	C
ATOM	858	CG ARG A 69	-68.332	27.445	-12.335	1.00	16.07	C
ATOM	861	CD ARG A 69	-69.378	26.578	-11.701	1.00	14.26	C
ATOM	864	NE ARG A 69	-68.915	25.205	-11.677	1.00	12.82	N
ATOM	866	CZ ARG A 69	-69.063	24.346	-12.676	1.00	12.42	C
ATOM	867	NH1 ARG A 69	-69.676	24.694	-13.792	1.00	12.27	N
ATOM	870	NH2 ARG A 69	-68.601	23.119	-12.553	1.00	13.23	N
ATOM	873	C ARG A 69	-66.010	29.409	-12.692	1.00	17.62	C
ATOM	874	O ARG A 69	-66.281	29.401	-13.892	1.00	17.10	O
ATOM	876	N PHE A 70	-65.811	30.529	-11.985	1.00	19.11	N
ATOM	877	CA PHE A 70	-65.936	31.891	-12.530	1.00	19.96	C
ATOM	879	CB PHE A 70	-67.024	32.666	-11.763	1.00	19.89	C
ATOM	882	CG PHE A 70	-68.365	32.010	-11.784	1.00	18.68	C
ATOM	883	CD1 PHE A 70	-69.158	32.082	-12.903	1.00	18.05	C
ATOM	885	CE1 PHE A 70	-70.376	31.452	-12.938	1.00	18.03	C
ATOM	887	CZ PHE A 70	-70.825	30.754	-11.851	1.00	17.26	C
ATOM	889	CE2 PHE A 70	-70.051	30.672	-10.736	1.00	17.60	C
ATOM	891	CD2 PHE A 70	-68.822	31.300	-10.700	1.00	17.86	C
ATOM	893	C PHE A 70	-64.627	32.661	-12.400	1.00	21.34	C
ATOM	894	O PHE A 70	-64.629	33.862	-12.171	1.00	21.25	O
ATOM	896	N GLU A 71	-63.501	31.977	-12.541	1.00	23.19	N
ATOM	897	CA GLU A 71	-62.211	32.591	-12.214	1.00	24.60	C
ATOM	899	CB GLU A 71	-61.127	31.524	-12.073	1.00	25.03	C
ATOM	902	CG GLU A 71	-59.717	32.059	-11.861	1.00	27.10	C
ATOM	905	CD GLU A 71	-58.712	30.949	-11.549	1.00	30.10	C
ATOM	906	OE1 GLU A 71	-59.084	30.007	-10.796	1.00	32.05	O
ATOM	907	OE2 GLU A 71	-57.558	31.027	-12.051	1.00	30.50	O

ATOM	908	C	GLU A 71	-61.809	33.644	-13.241	1.00	25.23	C
ATOM	909	O	GLU A 71	-61.362	34.717	-12.862	1.00	25.54	O
ATOM	911	N	SER A 72	-61.979	33.345	-14.530	1.00	25.91	N
ATOM	912	CA	SER A 72	-61.641	34.294	-15.588	1.00	26.37	C
ATOM	914	CB	SER A 72	-61.656	33.606	-16.941	1.00	26.32	C
ATOM	917	OG	SER A 72	-62.985	33.278	-17.291	1.00	27.28	O
ATOM	919	C	SER A 72	-62.609	35.476	-15.608	1.00	26.84	C
ATOM	920	O	SER A 72	-62.197	36.601	-15.852	1.00	27.02	O
ATOM	922	N	ASP A 73	-63.893	35.214	-15.361	1.00	27.45	N
ATOM	923	CA	ASP A 73	-64.895	36.277	-15.246	1.00	27.87	C
ATOM	925	CB	ASP A 73	-66.289	35.705	-14.968	1.00	27.99	C
ATOM	928	CG	ASP A 73	-66.842	34.913	-16.140	1.00	29.49	C
ATOM	929	OD1	ASP A 73	-67.215	35.530	-17.155	1.00	30.90	O
ATOM	930	OD2	ASP A 73	-66.916	33.667	-16.054	1.00	31.85	O
ATOM	931	C	ASP A 73	-64.518	37.221	-14.121	1.00	27.86	C
ATOM	932	O	ASP A 73	-64.598	38.435	-14.286	1.00	28.05	O
ATOM	934	N	ILE A 74	-64.118	36.652	-12.985	1.00	27.88	N
ATOM	935	CA	ILE A 74	-63.698	37.422	-11.821	1.00	27.98	C
ATOM	937	CB	ILE A 74	-63.348	36.521	-10.637	1.00	27.71	C
ATOM	939	CG1	ILE A 74	-64.607	36.007	-9.960	1.00	27.60	C
ATOM	942	CD1	ILE A 74	-64.355	34.846	-9.023	1.00	27.67	C
ATOM	946	CG2	ILE A 74	-62.551	37.272	-9.621	1.00	26.82	C
ATOM	950	C	ILE A 74	-62.472	38.252	-12.133	1.00	28.80	C
ATOM	951	O	ILE A 74	-62.475	39.453	-11.917	1.00	29.00	O
ATOM	953	N	ARG A 75	-61.415	37.616	-12.628	1.00	29.80	N
ATOM	954	CA	ARG A 75	-60.197	38.341	-12.968	1.00	30.79	C
ATOM	956	CB	ARG A 75	-59.286	37.507	-13.841	1.00	31.33	C
ATOM	959	CG	ARG A 75	-58.506	36.441	-13.115	1.00	34.25	C
ATOM	962	CD	ARG A 75	-57.286	36.009	-13.929	1.00	37.53	C
ATOM	965	NE	ARG A 75	-56.238	37.019	-13.799	1.00	40.89	N
ATOM	967	CZ	ARG A 75	-54.934	36.807	-13.967	1.00	44.24	C
ATOM	968	NH1	ARG A 75	-54.459	35.602	-14.294	1.00	45.02	N
ATOM	971	NH2	ARG A 75	-54.089	37.821	-13.796	1.00	45.88	N
ATOM	974	C	ARG A 75	-60.548	39.596	-13.727	1.00	30.95	C
ATOM	975	O	ARG A 75	-60.163	40.685	-13.336	1.00	31.27	O
ATOM	977	N	ARG A 76	-61.293	39.435	-14.815	1.00	31.31	N
ATOM	978	CA	ARG A 76	-61.712	40.567	-15.638	1.00	31.62	C
ATOM	980	CB	ARG A 76	-62.593	40.106	-16.794	1.00	31.97	C
ATOM	983	CG	ARG A 76	-61.833	39.392	-17.895	1.00	32.91	C
ATOM	986	CD	ARG A 76	-62.615	39.453	-19.205	1.00	34.43	C
ATOM	989	NE	ARG A 76	-63.925	38.803	-19.128	1.00	35.66	N
ATOM	991	CZ	ARG A 76	-64.115	37.480	-19.106	1.00	37.15	C
ATOM	992	NH1	ARG A 76	-63.085	36.632	-19.129	1.00	37.93	N
ATOM	995	NH2	ARG A 76	-65.347	36.992	-19.048	1.00	37.64	N
ATOM	998	C	ARG A 76	-62.453	41.629	-14.849	1.00	31.41	C
ATOM	999	O	ARG A 76	-62.144	42.798	-14.966	1.00	31.35	O
ATOM	1001	N	ALA A 77	-63.437	41.230	-14.058	1.00	31.53	N
ATOM	1002	CA	ALA A 77	-64.136	42.177	-13.211	1.00	31.73	C
ATOM	1004	CB	ALA A 77	-65.074	41.468	-12.275	1.00	31.69	C
ATOM	1008	C	ALA A 77	-63.107	42.950	-12.423	1.00	32.23	C
ATOM	1009	O	ALA A 77	-63.116	44.176	-12.411	1.00	32.42	O
ATOM	1011	N	LEU A 78	-62.198	42.223	-11.786	1.00	32.95	N
ATOM	1012	CA	LEU A 78	-61.168	42.838	-10.963	1.00	33.51	C
ATOM	1014	CB	LEU A 78	-60.319	41.775	-10.265	1.00	33.10	C
ATOM	1017	CG	LEU A 78	-60.959	41.077	-9.081	1.00	32.30	C
ATOM	1019	CD1	LEU A 78	-59.960	40.104	-8.472	1.00	30.80	C
ATOM	1023	CD2	LEU A 78	-61.438	42.117	-8.056	1.00	31.89	C
ATOM	1027	C	LEU A 78	-60.249	43.752	-11.753	1.00	34.70	C
ATOM	1028	O	LEU A 78	-59.790	44.753	-11.214	1.00	35.18	O
ATOM	1030	N	ASP A 79	-59.955	43.404	-13.007	1.00	35.91	N

ATOM	1031	CA	ASP	A	79	-59.054	44.205	-13.828	1.00	36.89	C
ATOM	1033	CB	ASP	A	79	-58.637	43.452	-15.083	1.00	37.06	C
ATOM	1036	CG	ASP	A	79	-57.518	44.150	-15.821	1.00	38.41	C
ATOM	1037	OD1	ASP	A	79	-57.800	45.080	-16.615	1.00	39.85	O
ATOM	1038	OD2	ASP	A	79	-56.347	43.775	-15.592	1.00	40.78	O
ATOM	1039	C	ASP	A	79	-59.721	45.505	-14.220	1.00	37.75	C
ATOM	1040	O	ASP	A	79	-59.105	46.568	-14.177	1.00	37.82	O
ATOM	1042	N	ARG	A	80	-60.984	45.401	-14.615	1.00	38.99	N
ATOM	1043	CA	ARG	A	80	-61.814	46.558	-14.923	1.00	39.88	C
ATOM	1045	CB	ARG	A	80	-63.191	46.106	-15.427	1.00	40.41	C
ATOM	1048	CG	ARG	A	80	-64.025	47.183	-16.122	1.00	42.77	C
ATOM	1051	CD	ARG	A	80	-65.347	46.624	-16.718	1.00	45.79	C
ATOM	1054	NE	ARG	A	80	-66.107	45.799	-15.765	1.00	48.55	N
ATOM	1056	CZ	ARG	A	80	-66.136	44.457	-15.742	1.00	50.37	C
ATOM	1057	NH1	ARG	A	80	-65.450	43.724	-16.636	1.00	49.78	N
ATOM	1060	NH2	ARG	A	80	-66.868	43.835	-14.805	1.00	50.77	N
ATOM	1063	C	ARG	A	80	-61.946	47.409	-13.668	1.00	39.80	C
ATOM	1064	O	ARG	A	80	-61.840	48.630	-13.742	1.00	40.37	O
ATOM	1066	N	PHE	A	81	-62.136	46.771	-12.516	1.00	39.62	N
ATOM	1067	CA	PHE	A	81	-62.216	47.505	-11.245	1.00	39.61	C
ATOM	1069	CB	PHE	A	81	-62.392	46.561	-10.053	1.00	39.74	C
ATOM	1072	CG	PHE	A	81	-62.282	47.247	-8.712	1.00	38.81	C
ATOM	1073	CD1	PHE	A	81	-63.238	48.153	-8.315	1.00	39.03	C
ATOM	1075	CE1	PHE	A	81	-63.153	48.787	-7.100	1.00	39.69	C
ATOM	1077	CZ	PHE	A	81	-62.095	48.515	-6.261	1.00	39.54	C
ATOM	1079	CE2	PHE	A	81	-61.134	47.608	-6.647	1.00	38.66	C
ATOM	1081	CD2	PHE	A	81	-61.228	46.986	-7.865	1.00	38.27	C
ATOM	1083	C	PHE	A	81	-61.006	48.378	-10.970	1.00	39.58	C
ATOM	1084	O	PHE	A	81	-61.165	49.535	-10.581	1.00	39.76	O
ATOM	1086	N	VAL	A	82	-59.810	47.817	-11.142	1.00	39.50	N
ATOM	1087	CA	VAL	A	82	-58.575	48.562	-10.892	1.00	39.47	C
ATOM	1089	CB	VAL	A	82	-57.314	47.659	-10.977	1.00	39.37	C
ATOM	1091	CG1	VAL	A	82	-56.108	48.426	-11.521	1.00	39.06	C
ATOM	1095	CG2	VAL	A	82	-57.009	47.079	-9.614	1.00	39.23	C
ATOM	1099	C	VAL	A	82	-58.457	49.755	-11.833	1.00	39.54	C
ATOM	1100	O	VAL	A	82	-58.222	50.884	-11.378	1.00	39.45	O
ATOM	1102	N	SER	A	83	-58.678	49.513	-13.126	1.00	39.65	N
ATOM	1103	CA	SER	A	83	-58.469	50.535	-14.157	1.00	39.82	C
ATOM	1105	CB	SER	A	83	-58.196	49.880	-15.526	1.00	39.84	C
ATOM	1108	OG	SER	A	83	-59.138	48.870	-15.829	1.00	40.00	O
ATOM	1110	C	SER	A	83	-59.617	51.561	-14.219	1.00	39.62	C
ATOM	1111	O	SER	A	83	-60.149	51.853	-15.286	1.00	39.59	O
ATOM	1113	N	SER	A	84	-59.957	52.108	-13.052	1.00	39.44	N
ATOM	1114	CA	SER	A	84	-60.926	53.197	-12.910	1.00	39.13	C
ATOM	1116	CB	SER	A	84	-62.238	52.834	-13.618	1.00	39.03	C
ATOM	1119	OG	SER	A	84	-62.620	51.505	-13.325	1.00	37.97	O
ATOM	1121	C	SER	A	84	-61.243	53.580	-11.446	1.00	39.17	C
ATOM	1122	O	SER	A	84	-62.189	54.341	-11.230	1.00	39.58	O
ATOM	1124	N	GLY	A	85	-60.482	53.077	-10.457	1.00	38.74	N
ATOM	1125	CA	GLY	A	85	-60.795	53.282	-9.024	1.00	38.26	C
ATOM	1128	C	GLY	A	85	-61.724	52.238	-8.413	1.00	37.79	C
ATOM	1129	O	GLY	A	85	-61.529	51.799	-7.273	1.00	36.84	O
ATOM	1131	N	SER	A	94	-63.557	54.187	1.994	1.00	28.27	N
ATOM	1132	CA	SER	A	94	-64.069	52.892	2.476	1.00	28.03	C
ATOM	1134	CB	SER	A	94	-64.856	52.225	1.359	1.00	27.84	C
ATOM	1137	OG	SER	A	94	-65.234	50.921	1.743	1.00	28.20	O
ATOM	1139	C	SER	A	94	-62.992	51.910	2.934	1.00	27.83	C
ATOM	1140	O	SER	A	94	-62.339	51.325	2.090	1.00	28.15	O
ATOM	1142	N	LEU	A	95	-62.815	51.705	4.245	1.00	27.70	N
ATOM	1143	CA	LEU	A	95	-61.813	50.728	4.751	1.00	27.59	C

ATOM	1145	CB	LEU	A	95	-61.692	50.747	6.282	1.00	27.49	C
ATOM	1148	CG	LEU	A	95	-60.871	49.596	6.909	1.00	27.47	C
ATOM	1150	CD1	LEU	A	95	-59.400	49.900	6.902	1.00	27.70	C
ATOM	1154	CD2	LEU	A	95	-61.292	49.274	8.329	1.00	27.50	C
ATOM	1158	C	LEU	A	95	-62.173	49.314	4.318	1.00	27.71	C
ATOM	1159	O	LEU	A	95	-61.303	48.537	3.914	1.00	27.73	O
ATOM	1161	N	HIS	A	96	-63.461	48.989	4.437	1.00	27.74	N
ATOM	1162	CA	HIS	A	96	-64.003	47.715	3.977	1.00	27.54	C
ATOM	1164	CB	HIS	A	96	-65.497	47.680	4.230	1.00	27.57	C
ATOM	1167	CG	HIS	A	96	-66.161	46.482	3.654	1.00	28.65	C
ATOM	1168	ND1	HIS	A	96	-65.790	45.201	3.988	1.00	30.59	N
ATOM	1170	CE1	HIS	A	96	-66.540	44.341	3.325	1.00	31.88	C
ATOM	1172	NE2	HIS	A	96	-67.384	45.021	2.570	1.00	32.06	N
ATOM	1174	CD2	HIS	A	96	-67.166	46.364	2.758	1.00	30.47	C
ATOM	1176	C	HIS	A	96	-63.705	47.450	2.492	1.00	27.09	C
ATOM	1177	O	HIS	A	96	-63.128	46.425	2.150	1.00	27.18	O
ATOM	1179	N	GLY	A	97	-64.079	48.386	1.625	1.00	26.55	N
ATOM	1180	CA	GLY	A	97	-63.760	48.300	.205	1.00	26.26	C
ATOM	1183	C	GLY	A	97	-62.277	48.112	-.088	1.00	26.23	C
ATOM	1184	O	GLY	A	97	-61.918	47.334	-.974	1.00	26.39	O
ATOM	1186	N	THR	A	98	-61.415	48.815	.653	1.00	25.83	N
ATOM	1187	CA	THR	A	98	-59.967	48.755	.435	1.00	25.36	C
ATOM	1189	CB	THR	A	98	-59.204	49.847	1.215	1.00	25.07	C
ATOM	1191	OG1	THR	A	98	-59.796	51.123	.980	1.00	23.62	O
ATOM	1193	CG2	THR	A	98	-57.755	49.905	.780	1.00	24.88	C
ATOM	1197	C	THR	A	98	-59.416	47.389	.845	1.00	25.69	C
ATOM	1198	O	THR	A	98	-58.622	46.802	.111	1.00	25.66	O
ATOM	1200	N	ALA	A	99	-59.847	46.890	2.008	1.00	25.95	N
ATOM	1201	CA	ALA	A	99	-59.468	45.540	2.495	1.00	25.93	C
ATOM	1203	CB	ALA	A	99	-59.945	45.330	3.929	1.00	25.82	C
ATOM	1207	C	ALA	A	99	-59.984	44.404	1.595	1.00	25.64	C
ATOM	1208	O	ALA	A	99	-59.204	43.599	1.102	1.00	25.38	O
ATOM	1210	N	LEU	A	100	-61.290	44.348	1.371	1.00	25.46	N
ATOM	1211	CA	LEU	A	100	-61.856	43.335	.474	1.00	25.55	C
ATOM	1213	CB	LEU	A	100	-63.361	43.540	.326	1.00	25.37	C
ATOM	1216	CG	LEU	A	100	-64.118	42.464	-.443	1.00	25.41	C
ATOM	1218	CD1	LEU	A	100	-64.038	41.115	.281	1.00	26.01	C
ATOM	1222	CD2	LEU	A	100	-65.555	42.893	-.639	1.00	25.19	C
ATOM	1226	C	LEU	A	100	-61.196	43.353	-.919	1.00	25.70	C
ATOM	1227	O	LEU	A	100	-60.819	42.301	-1.469	1.00	25.36	O
ATOM	1229	N	SER	A	101	-61.064	44.550	-1.489	1.00	25.77	N
ATOM	1230	CA	SER	A	101	-60.434	44.688	-2.805	1.00	25.67	C
ATOM	1232	CB	SER	A	101	-60.629	46.098	-3.375	1.00	25.83	C
ATOM	1235	OG	SER	A	101	-59.840	47.064	-2.688	1.00	26.98	O
ATOM	1237	C	SER	A	101	-58.947	44.336	-2.756	1.00	25.01	C
ATOM	1238	O	SER	A	101	-58.427	43.705	-3.672	1.00	24.94	O
ATOM	1240	N	PHE	A	102	-58.270	44.737	-1.687	1.00	24.30	N
ATOM	1241	CA	PHE	A	102	-56.852	44.436	-1.560	1.00	23.99	C
ATOM	1243	CB	PHE	A	102	-56.295	44.963	-.229	1.00	23.84	C
ATOM	1246	CG	PHE	A	102	-54.860	44.606	.027	1.00	23.24	C
ATOM	1247	CD1	PHE	A	102	-53.849	45.504	-.259	1.00	23.56	C
ATOM	1249	CE1	PHE	A	102	-52.507	45.182	-.007	1.00	24.50	C
ATOM	1251	CZ	PHE	A	102	-52.175	43.941	.540	1.00	24.27	C
ATOM	1253	CE2	PHE	A	102	-53.180	43.042	.834	1.00	23.95	C
ATOM	1255	CD2	PHE	A	102	-54.520	43.376	.576	1.00	23.45	C
ATOM	1257	C	PHE	A	102	-56.673	42.928	-1.679	1.00	23.65	C
ATOM	1258	O	PHE	A	102	-55.890	42.442	-2.501	1.00	23.55	O
ATOM	1260	N	ARG	A	103	-57.443	42.202	-.880	1.00	23.20	N
ATOM	1261	CA	ARG	A	103	-57.298	40.765	-.764	1.00	22.84	C
ATOM	1263	CB	ARG	A	103	-58.250	40.231	.291	1.00	23.01	C

ATOM	1266	CG	ARG A 103	-58.054	38.774	.597	1.00	23.93	C
ATOM	1269	CD	ARG A 103	-58.703	38.376	1.937	1.00	23.84	C
ATOM	1272	NE	ARG A 103	-60.150	38.538	1.919	1.00	22.36	N
ATOM	1274	CZ	ARG A 103	-60.982	37.763	1.240	1.00	20.62	C
ATOM	1275	NH1	ARG A 103	-60.512	36.773	.490	1.00	20.19	N
ATOM	1278	NH2	ARG A 103	-62.288	38.003	1.293	1.00	20.32	N
ATOM	1281	C	ARG A 103	-57.602	40.117	-2.082	1.00	22.34	C
ATOM	1282	O	ARG A 103	-56.826	39.303	-2.576	1.00	22.37	O
ATOM	1284	N	LEU A 104	-58.729	40.496	-2.667	1.00	21.98	N
ATOM	1285	CA	LEU A 104	-59.129	39.926	-3.951	1.00	21.60	C
ATOM	1287	CB	LEU A 104	-60.503	40.442	-4.378	1.00	21.23	C
ATOM	1290	CG	LEU A 104	-61.590	39.832	-3.522	1.00	20.30	C
ATOM	1292	CD1	LEU A 104	-62.949	40.401	-3.856	1.00	19.53	C
ATOM	1296	CD2	LEU A 104	-61.554	38.344	-3.725	1.00	19.71	C
ATOM	1300	C	LEU A 104	-58.101	40.198	-5.040	1.00	21.48	C
ATOM	1301	O	LEU A 104	-57.743	39.294	-5.780	1.00	21.62	O
ATOM	1303	N	LEU A 105	-57.623	41.432	-5.139	1.00	21.14	N
ATOM	1304	CA	LEU A 105	-56.655	41.747	-6.173	1.00	21.05	C
ATOM	1306	CB	LEU A 105	-56.352	43.248	-6.205	1.00	20.94	C
ATOM	1309	CG	LEU A 105	-57.465	44.165	-6.732	1.00	20.06	C
ATOM	1311	CD1	LEU A 105	-57.060	45.617	-6.580	1.00	17.92	C
ATOM	1315	CD2	LEU A 105	-57.804	43.857	-8.174	1.00	19.18	C
ATOM	1319	C	LEU A 105	-55.372	40.908	-5.997	1.00	21.32	C
ATOM	1320	O	LEU A 105	-54.840	40.357	-6.976	1.00	21.37	O
ATOM	1322	N	ARG A 106	-54.893	40.777	-4.763	1.00	21.25	N
ATOM	1323	CA	ARG A 106	-53.678	40.003	-4.528	1.00	21.33	C
ATOM	1325	CB	ARG A 106	-53.151	40.211	-3.117	1.00	21.53	C
ATOM	1328	CG	ARG A 106	-51.772	39.588	-2.917	1.00	22.60	C
ATOM	1331	CD	ARG A 106	-51.098	40.091	-1.666	1.00	23.60	C
ATOM	1334	NE	ARG A 106	-50.374	41.330	-1.897	1.00	24.60	N
ATOM	1336	CZ	ARG A 106	-49.606	41.916	-.989	1.00	26.54	C
ATOM	1337	NH1	ARG A 106	-49.469	41.380	.225	1.00	26.98	N
ATOM	1340	NH2	ARG A 106	-48.967	43.041	-1.295	1.00	27.19	N
ATOM	1343	C	ARG A 106	-53.888	38.517	-4.767	1.00	21.19	C
ATOM	1344	O	ARG A 106	-52.998	37.821	-5.278	1.00	21.18	O
ATOM	1346	N	GLN A 107	-55.061	38.028	-4.378	1.00	21.09	N
ATOM	1347	CA	GLN A 107	-55.425	36.635	-4.618	1.00	20.76	C
ATOM	1349	CB	GLN A 107	-56.861	36.378	-4.168	1.00	20.51	C
ATOM	1352	CG	GLN A 107	-57.329	34.941	-4.327	1.00	19.68	C
ATOM	1355	CD	GLN A 107	-58.790	34.772	-4.010	1.00	18.26	C
ATOM	1356	OE1	GLN A 107	-59.395	35.591	-3.334	1.00	17.44	O
ATOM	1357	NE2	GLN A 107	-59.366	33.698	-4.497	1.00	18.48	N
ATOM	1360	C	GLN A 107	-55.304	36.306	-6.094	1.00	20.93	C
ATOM	1361	O	GLN A 107	-54.917	35.212	-6.447	1.00	21.45	O
ATOM	1363	N	HIS A 108	-55.642	37.260	-6.951	1.00	21.09	N
ATOM	1364	CA	HIS A 108	-55.686	37.030	-8.379	1.00	21.28	C
ATOM	1366	CB	HIS A 108	-57.024	37.541	-8.913	1.00	21.02	C
ATOM	1369	CG	HIS A 108	-58.182	36.706	-8.478	1.00	19.70	C
ATOM	1370	ND1	HIS A 108	-58.637	35.632	-9.210	1.00	19.91	N
ATOM	1372	CE1	HIS A 108	-59.641	35.060	-8.571	1.00	19.27	C
ATOM	1374	NE2	HIS A 108	-59.851	35.720	-7.449	1.00	18.10	N
ATOM	1376	CD2	HIS A 108	-58.950	36.753	-7.367	1.00	18.98	C
ATOM	1378	C	HIS A 108	-54.456	37.621	-9.108	1.00	22.21	C
ATOM	1379	O	HIS A 108	-54.505	37.980	-10.304	1.00	21.80	O
ATOM	1381	N	GLY A 109	-53.345	37.686	-8.375	1.00	23.19	N
ATOM	1382	CA	GLY A 109	-52.055	38.007	-8.956	1.00	24.17	C
ATOM	1385	C	GLY A 109	-51.813	39.470	-9.275	1.00	25.25	C
ATOM	1386	O	GLY A 109	-50.767	39.802	-9.845	1.00	25.26	O
ATOM	1388	N	PHE A 110	-52.755	40.349	-8.917	1.00	26.42	N
ATOM	1389	CA	PHE A 110	-52.575	41.783	-9.154	1.00	27.33	C

ATOM	1391	CB	PHE A 110	-53.895	42.546	-9.067	1.00	27.36	C
ATOM	1394	CG	PHE A 110	-54.838	42.250	-10.193	1.00	27.83	C
ATOM	1395	CD1	PHE A 110	-54.633	42.801	-11.439	1.00	28.69	C
ATOM	1397	CE1	PHE A 110	-55.490	42.529	-12.488	1.00	28.91	C
ATOM	1399	CZ	PHE A 110	-56.568	41.702	-12.295	1.00	28.85	C
ATOM	1401	CE2	PHE A 110	-56.787	41.146	-11.060	1.00	28.67	C
ATOM	1403	CD2	PHE A 110	-55.923	41.418	-10.013	1.00	28.39	C
ATOM	1405	C	PHE A 110	-51.582	42.369	-8.163	1.00	28.25	C
ATOM	1406	O	PHE A 110	-51.449	41.890	-7.024	1.00	28.71	O
ATOM	1408	N	GLU A 111	-50.886	43.411	-8.607	1.00	29.05	N
ATOM	1409	CA	GLU A 111	-49.951	44.142	-7.760	1.00	29.66	C
ATOM	1411	CB	GLU A 111	-48.902	44.803	-8.648	1.00	30.37	C
ATOM	1414	CG	GLU A 111	-47.750	45.509	-7.947	1.00	33.30	C
ATOM	1417	CD	GLU A 111	-46.916	46.342	-8.946	1.00	38.00	C
ATOM	1418	OE1	GLU A 111	-46.312	45.750	-9.885	1.00	39.60	O
ATOM	1419	OE2	GLU A 111	-46.886	47.594	-8.802	1.00	40.87	O
ATOM	1420	C	GLU A 111	-50.724	45.184	-6.935	1.00	28.98	C
ATOM	1421	O	GLU A 111	-51.293	46.123	-7.495	1.00	28.45	O
ATOM	1423	N	VAL A 112	-50.787	44.961	-5.618	1.00	28.47	N
ATOM	1424	CA	VAL A 112	-51.284	45.950	-4.655	1.00	28.09	C
ATOM	1426	CB	VAL A 112	-52.564	45.501	-3.908	1.00	27.76	C
ATOM	1428	CG1	VAL A 112	-53.745	45.466	-4.836	1.00	27.38	C
ATOM	1432	CG2	VAL A 112	-52.360	44.167	-3.252	1.00	27.80	C
ATOM	1436	C	VAL A 112	-50.193	46.199	-3.631	1.00	28.15	C
ATOM	1437	O	VAL A 112	-49.339	45.354	-3.436	1.00	27.90	O
ATOM	1439	N	SER A 113	-50.241	47.354	-2.974	1.00	28.54	N
ATOM	1440	CA	SER A 113	-49.194	47.782	-2.047	1.00	28.93	C
ATOM	1442	CB	SER A 113	-48.563	49.064	-2.564	1.00	29.12	C
ATOM	1445	OG	SER A 113	-47.692	49.644	-1.619	1.00	29.57	O
ATOM	1447	C	SER A 113	-49.780	48.041	-.685	1.00	29.13	C
ATOM	1448	O	SER A 113	-50.937	48.426	-.581	1.00	28.90	O
ATOM	1450	N	GLN A 114	-48.985	47.857	.367	1.00	29.84	N
ATOM	1451	CA	GLN A 114	-49.524	48.019	1.736	1.00	30.40	C
ATOM	1453	CB	GLN A 114	-48.568	47.507	2.830	1.00	30.34	C
ATOM	1456	CG	GLN A 114	-47.199	48.159	2.844	1.00	30.61	C
ATOM	1459	CD	GLN A 114	-46.456	47.928	4.147	1.00	30.77	C
ATOM	1460	OE1	GLN A 114	-46.876	47.129	4.982	1.00	31.92	O
ATOM	1461	NE2	GLN A 114	-45.341	48.624	4.324	1.00	30.02	N
ATOM	1464	C	GLN A 114	-49.939	49.440	2.068	1.00	30.62	C
ATOM	1465	O	GLN A 114	-50.573	49.644	3.094	1.00	30.73	O
ATOM	1467	N	GLU A 115	-49.592	50.407	1.216	1.00	30.94	N
ATOM	1468	CA	GLU A 115	-49.963	51.794	1.458	1.00	31.45	C
ATOM	1470	CB	GLU A 115	-48.881	52.777	.975	1.00	31.86	C
ATOM	1473	CG	GLU A 115	-48.489	52.723	-.509	1.00	33.48	C
ATOM	1476	CD	GLU A 115	-47.006	53.100	-.759	1.00	36.13	C
ATOM	1477	OE1	GLU A 115	-46.106	52.668	.022	1.00	38.72	O
ATOM	1478	OE2	GLU A 115	-46.741	53.811	-1.754	1.00	35.70	O
ATOM	1479	C	GLU A 115	-51.338	52.125	.895	1.00	31.45	C
ATOM	1480	O	GLU A 115	-51.716	53.283	.808	1.00	31.80	O
ATOM	1482	N	ALA A 116	-52.108	51.105	.549	1.00	31.58	N
ATOM	1483	CA	ALA A 116	-53.524	51.286	.284	1.00	31.67	C
ATOM	1485	CB	ALA A 116	-54.071	50.092	-.461	1.00	31.62	C
ATOM	1489	C	ALA A 116	-54.273	51.473	1.596	1.00	31.88	C
ATOM	1490	O	ALA A 116	-55.428	51.874	1.604	1.00	31.82	O
ATOM	1492	N	PHE A 117	-53.608	51.161	2.702	1.00	32.39	N
ATOM	1493	CA	PHE A 117	-54.193	51.271	4.030	1.00	32.87	C
ATOM	1495	CB	PHE A 117	-53.884	50.003	4.856	1.00	32.79	C
ATOM	1498	CG	PHE A 117	-54.539	48.758	4.323	1.00	31.03	C
ATOM	1499	CD1	PHE A 117	-53.782	47.749	3.760	1.00	28.79	C
ATOM	1501	CE1	PHE A 117	-54.379	46.620	3.265	1.00	28.36	C

ATOM	1503	CZ	PHE A 117	-55.754	46.489	3.316	1.00	29.50	C
ATOM	1505	CE2	PHE A 117	-56.527	47.493	3.872	1.00	29.91	C
ATOM	1507	CD2	PHE A 117	-55.917	48.613	4.375	1.00	30.09	C
ATOM	1509	C	PHE A 117	-53.683	52.493	4.773	1.00	33.71	C
ATOM	1510	O	PHE A 117	-54.012	52.678	5.937	1.00	33.61	O
ATOM	1512	N	SER A 118	-52.886	53.325	4.107	1.00	35.09	N
ATOM	1513	CA	SER A 118	-52.186	54.429	4.779	1.00	36.12	C
ATOM	1515	CB	SER A 118	-51.073	54.990	3.895	1.00	36.09	C
ATOM	1518	OG	SER A 118	-51.607	55.747	2.820	1.00	35.96	O
ATOM	1520	C	SER A 118	-53.134	55.552	5.184	1.00	37.17	C
ATOM	1521	O	SER A 118	-52.873	56.255	6.166	1.00	37.22	O
ATOM	1523	N	GLY A 119	-54.231	55.705	4.436	1.00	38.40	N
ATOM	1524	CA	GLY A 119	-55.227	56.735	4.708	1.00	39.47	C
ATOM	1527	C	GLY A 119	-56.255	56.398	5.778	1.00	40.62	C
ATOM	1528	O	GLY A 119	-57.259	57.086	5.877	1.00	40.74	O
ATOM	1530	N	PHE A 120	-56.024	55.344	6.563	1.00	42.13	N
ATOM	1531	CA	PHE A 120	-56.906	54.971	7.684	1.00	43.27	C
ATOM	1533	CB	PHE A 120	-57.510	53.586	7.455	1.00	43.21	C
ATOM	1536	CG	PHE A 120	-58.176	53.439	6.134	1.00	42.36	C
ATOM	1537	CD1	PHE A 120	-59.418	54.006	5.911	1.00	41.23	C
ATOM	1539	CE1	PHE A 120	-60.035	53.880	4.691	1.00	40.87	C
ATOM	1541	CZ	PHE A 120	-59.414	53.182	3.675	1.00	41.03	C
ATOM	1543	CE2	PHE A 120	-58.167	52.618	3.883	1.00	40.92	C
ATOM	1545	CD2	PHE A 120	-57.556	52.746	5.106	1.00	41.30	C
ATOM	1547	C	PHE A 120	-56.179	54.952	9.019	1.00	44.71	C
ATOM	1548	O	PHE A 120	-56.732	54.489	10.025	1.00	44.80	O
ATOM	1550	N	LYS A 121	-54.939	55.430	9.021	1.00	46.38	N
ATOM	1551	CA	LYS A 121	-54.153	55.519	10.235	1.00	47.81	C
ATOM	1553	CB	LYS A 121	-52.722	55.050	9.957	1.00	47.86	C
ATOM	1556	CG	LYS A 121	-52.648	53.537	9.648	1.00	48.24	C
ATOM	1559	CD	LYS A 121	-51.417	53.119	8.823	1.00	49.41	C
ATOM	1562	CE	LYS A 121	-50.095	53.142	9.621	1.00	49.88	C
ATOM	1565	NZ	LYS A 121	-49.944	52.005	10.584	1.00	49.89	N
ATOM	1569	C	LYS A 121	-54.257	56.964	10.735	1.00	49.00	C
ATOM	1570	O	LYS A 121	-54.712	57.839	10.000	1.00	49.11	O
ATOM	1572	N	ASP A 122	-53.899	57.206	11.995	1.00	50.51	N
ATOM	1573	CA	ASP A 122	-54.078	58.537	12.596	1.00	51.46	C
ATOM	1575	CB	ASP A 122	-54.604	58.443	14.050	1.00	51.29	C
ATOM	1578	CG	ASP A 122	-53.570	57.919	15.036	1.00	50.73	C
ATOM	1579	OD1	ASP A 122	-52.353	58.025	14.777	1.00	49.90	O
ATOM	1580	OD2	ASP A 122	-53.988	57.398	16.088	1.00	49.66	O
ATOM	1581	C	ASP A 122	-52.792	59.358	12.502	1.00	52.60	C
ATOM	1582	O	ASP A 122	-51.830	58.953	11.835	1.00	52.65	O
ATOM	1584	N	GLN A 123	-52.802	60.512	13.172	1.00	53.88	N
ATOM	1585	CA	GLN A 123	-51.680	61.456	13.195	1.00	54.66	C
ATOM	1587	CB	GLN A 123	-51.993	62.594	14.178	1.00	55.01	C
ATOM	1590	CG	GLN A 123	-53.249	63.442	13.831	1.00	56.15	C
ATOM	1593	CD	GLN A 123	-52.943	64.736	13.064	1.00	57.64	C
ATOM	1594	OE1	GLN A 123	-51.779	65.102	12.859	1.00	59.10	O
ATOM	1595	NE2	GLN A 123	-54.001	65.438	12.652	1.00	57.63	N
ATOM	1598	C	GLN A 123	-50.365	60.772	13.591	1.00	54.87	C
ATOM	1599	O	GLN A 123	-49.309	61.078	13.034	1.00	54.67	O
ATOM	1601	N	ASN A 124	-50.458	59.834	14.538	1.00	55.21	N
ATOM	1602	CA	ASN A 124	-49.300	59.137	15.109	1.00	55.39	C
ATOM	1604	CB	ASN A 124	-49.421	59.122	16.637	1.00	55.52	C
ATOM	1607	CG	ASN A 124	-49.833	60.478	17.202	1.00	55.98	C
ATOM	1608	OD1	ASN A 124	-49.093	61.460	17.091	1.00	56.93	O
ATOM	1609	ND2	ASN A 124	-51.026	60.540	17.797	1.00	55.99	N
ATOM	1612	C	ASN A 124	-49.116	57.705	14.589	1.00	55.27	C
ATOM	1613	O	ASN A 124	-48.530	56.864	15.271	1.00	55.18	O

ATOM	1615	N	GLY A 125	-49.626	57.430	13.391	1.00	55.18	N
ATOM	1616	CA	GLY A 125	-49.311	56.190	12.672	1.00	55.08	C
ATOM	1619	C	GLY A 125	-49.994	54.897	13.119	1.00	54.90	C
ATOM	1620	O	GLY A 125	-49.543	53.806	12.750	1.00	55.24	O
ATOM	1622	N	ASN A 126	-51.072	55.008	13.900	1.00	54.23	N
ATOM	1623	CA	ASN A 126	-51.871	53.853	14.329	1.00	53.39	C
ATOM	1625	CB	ASN A 126	-51.927	53.782	15.852	1.00	53.27	C
ATOM	1628	CG	ASN A 126	-50.596	53.425	16.466	1.00	52.81	C
ATOM	1629	OD1	ASN A 126	-49.636	53.113	15.764	1.00	52.29	O
ATOM	1630	ND2	ASN A 126	-50.532	53.460	17.789	1.00	52.53	N
ATOM	1633	C	ASN A 126	-53.274	53.978	13.771	1.00	52.82	C
ATOM	1634	O	ASN A 126	-53.724	55.081	13.494	1.00	52.69	O
ATOM	1636	N	PHE A 127	-53.968	52.857	13.607	1.00	52.17	N
ATOM	1637	CA	PHE A 127	-55.287	52.875	12.962	1.00	51.56	C
ATOM	1639	CB	PHE A 127	-55.805	51.451	12.721	1.00	51.35	C
ATOM	1642	CG	PHE A 127	-55.119	50.730	11.587	1.00	50.06	C
ATOM	1643	CD1	PHE A 127	-54.068	49.856	11.831	1.00	48.90	C
ATOM	1645	CE1	PHE A 127	-53.437	49.187	10.789	1.00	47.59	C
ATOM	1647	CZ	PHE A 127	-53.857	49.383	9.495	1.00	47.40	C
ATOM	1649	CE2	PHE A 127	-54.906	50.244	9.234	1.00	47.96	C
ATOM	1651	CD2	PHE A 127	-55.533	50.915	10.276	1.00	48.85	C
ATOM	1653	C	PHE A 127	-56.303	53.690	13.776	1.00	51.35	C
ATOM	1654	O	PHE A 127	-56.347	53.594	15.002	1.00	51.09	O
ATOM	1656	N	LEU A 128	-57.100	54.500	13.083	1.00	51.15	N
ATOM	1657	CA	LEU A 128	-58.143	55.280	13.727	1.00	51.09	C
ATOM	1659	CB	LEU A 128	-59.039	55.979	12.694	1.00	51.17	C
ATOM	1662	CG	LEU A 128	-58.463	57.046	11.746	1.00	51.31	C
ATOM	1664	CD1	LEU A 128	-59.571	57.628	10.855	1.00	50.76	C
ATOM	1668	CD2	LEU A 128	-57.762	58.157	12.508	1.00	50.98	C
ATOM	1672	C	LEU A 128	-58.996	54.353	14.581	1.00	51.06	C
ATOM	1673	O	LEU A 128	-59.573	53.393	14.078	1.00	50.94	O
ATOM	1675	N	GLU A 129	-59.053	54.637	15.879	1.00	51.12	N
ATOM	1676	CA	GLU A 129	-59.916	53.913	16.808	1.00	50.99	C
ATOM	1678	CB	GLU A 129	-60.019	54.665	18.132	1.00	51.14	C
ATOM	1681	CG	GLU A 129	-58.850	54.431	19.049	1.00	52.04	C
ATOM	1684	CD	GLU A 129	-58.903	53.069	19.697	1.00	53.22	C
ATOM	1685	OE1	GLU A 129	-58.040	52.219	19.376	1.00	54.19	O
ATOM	1686	OE2	GLU A 129	-59.823	52.848	20.516	1.00	53.51	O
ATOM	1687	C	GLU A 129	-61.314	53.712	16.260	1.00	50.66	C
ATOM	1688	O	GLU A 129	-61.791	52.587	16.200	1.00	50.83	O
ATOM	1690	N	ASN A 130	-61.960	54.799	15.839	1.00	50.28	N
ATOM	1691	CA	ASN A 130	-63.389	54.765	15.486	1.00	50.03	C
ATOM	1693	CB	ASN A 130	-63.910	56.182	15.193	1.00	50.09	C
ATOM	1696	CG	ASN A 130	-63.286	56.805	13.960	1.00	50.37	C
ATOM	1697	OD1	ASN A 130	-62.504	56.177	13.249	1.00	51.44	O
ATOM	1698	ND2	ASN A 130	-63.634	58.056	13.701	1.00	50.27	N
ATOM	1701	C	ASN A 130	-63.804	53.782	14.369	1.00	49.56	C
ATOM	1702	O	ASN A 130	-64.993	53.581	14.131	1.00	49.44	O
ATOM	1704	N	LEU A 131	-62.828	53.175	13.699	1.00	49.21	N
ATOM	1705	CA	LEU A 131	-63.087	52.103	12.729	1.00	48.78	C
ATOM	1707	CB	LEU A 131	-61.846	51.832	11.875	1.00	48.64	C
ATOM	1710	CG	LEU A 131	-61.445	52.939	10.898	1.00	48.26	C
ATOM	1712	CD1	LEU A 131	-60.029	52.709	10.375	1.00	47.76	C
ATOM	1716	CD2	LEU A 131	-62.441	53.039	9.754	1.00	47.32	C
ATOM	1720	C	LEU A 131	-63.516	50.798	13.394	1.00	48.58	C
ATOM	1721	O	LEU A 131	-64.093	49.939	12.731	1.00	48.64	O
ATOM	1723	N	LYS A 132	-63.225	50.645	14.689	1.00	48.23	N
ATOM	1724	CA	LYS A 132	-63.650	49.474	15.465	1.00	47.80	C
ATOM	1726	CB	LYS A 132	-63.156	49.587	16.913	1.00	47.71	C
ATOM	1729	CG	LYS A 132	-63.930	50.615	17.728	1.00	47.80	C

ATOM	1732	CD	LYS A 132	-63.430	50.779	19.153	1.00	47.75	C
ATOM	1735	CE	LYS A 132	-64.438	51.577	19.982	1.00	47.22	C
ATOM	1738	NZ	LYS A 132	-63.769	52.364	21.033	1.00	46.65	N
ATOM	1742	C	LYS A 132	-65.182	49.308	15.446	1.00	47.56	C
ATOM	1743	O	LYS A 132	-65.696	48.204	15.616	1.00	47.76	O
ATOM	1745	N	GLU A 133	-65.901	50.413	15.249	1.00	47.18	N
ATOM	1746	CA	GLU A 133	-67.371	50.412	15.191	1.00	46.79	C
ATOM	1748	CB	GLU A 133	-67.898	51.850	15.337	1.00	46.95	C
ATOM	1751	CG	GLU A 133	-67.840	52.376	16.776	1.00	47.56	C
ATOM	1754	CD	GLU A 133	-67.613	53.886	16.863	1.00	48.41	C
ATOM	1755	OE1	GLU A 133	-68.247	54.648	16.093	1.00	47.67	O
ATOM	1756	OE2	GLU A 133	-66.798	54.304	17.721	1.00	48.92	O
ATOM	1757	C	GLU A 133	-67.949	49.770	13.918	1.00	46.10	C
ATOM	1758	O	GLU A 133	-69.140	49.446	13.874	1.00	46.11	O
ATOM	1760	N	ASP A 134	-67.103	49.601	12.896	1.00	45.11	N
ATOM	1761	CA	ASP A 134	-67.477	48.995	11.608	1.00	44.06	C
ATOM	1763	CB	ASP A 134	-66.934	49.868	10.474	1.00	43.97	C
ATOM	1766	CG	ASP A 134	-67.380	49.411	9.104	1.00	44.08	C
ATOM	1767	OD1	ASP A 134	-68.072	48.373	8.972	1.00	43.42	O
ATOM	1768	OD2	ASP A 134	-67.020	50.116	8.141	1.00	44.84	O
ATOM	1769	C	ASP A 134	-66.910	47.570	11.526	1.00	43.12	C
ATOM	1770	O	ASP A 134	-65.788	47.357	11.080	1.00	43.19	O
ATOM	1772	N	ILE A 135	-67.708	46.594	11.944	1.00	41.85	N
ATOM	1773	CA	ILE A 135	-67.216	45.247	12.216	1.00	40.75	C
ATOM	1775	CB	ILE A 135	-68.183	44.507	13.170	1.00	40.80	C
ATOM	1777	CG1	ILE A 135	-68.227	45.226	14.524	1.00	41.21	C
ATOM	1780	CD1	ILE A 135	-69.353	46.283	14.633	1.00	42.86	C
ATOM	1784	CG2	ILE A 135	-67.776	43.063	13.373	1.00	40.88	C
ATOM	1788	C	ILE A 135	-66.980	44.457	10.933	1.00	39.67	C
ATOM	1789	O	ILE A 135	-66.166	43.537	10.901	1.00	39.37	O
ATOM	1791	N	LYS A 136	-67.690	44.828	9.877	1.00	38.62	N
ATOM	1792	CA	LYS A 136	-67.465	44.256	8.554	1.00	37.88	C
ATOM	1794	CB	LYS A 136	-68.525	44.749	7.561	1.00	38.33	C
ATOM	1797	CG	LYS A 136	-69.955	44.342	7.879	1.00	40.17	C
ATOM	1800	CD	LYS A 136	-70.486	43.262	6.911	1.00	42.97	C
ATOM	1803	CE	LYS A 136	-71.624	42.389	7.536	1.00	44.00	C
ATOM	1806	NZ	LYS A 136	-72.950	43.086	7.726	1.00	43.82	N
ATOM	1810	C	LYS A 136	-66.097	44.687	8.047	1.00	36.38	C
ATOM	1811	O	LYS A 136	-65.399	43.914	7.388	1.00	36.24	O
ATOM	1813	N	ALA A 137	-65.728	45.931	8.341	1.00	34.51	N
ATOM	1814	CA	ALA A 137	-64.463	46.474	7.876	1.00	33.41	C
ATOM	1816	CB	ALA A 137	-64.432	47.981	8.058	1.00	33.43	C
ATOM	1820	C	ALA A 137	-63.291	45.827	8.600	1.00	32.30	C
ATOM	1821	O	ALA A 137	-62.288	45.481	7.995	1.00	31.89	O
ATOM	1823	N	ILE A 138	-63.437	45.662	9.905	1.00	31.41	N
ATOM	1824	CA	ILE A 138	-62.395	45.074	10.737	1.00	30.60	C
ATOM	1826	CB	ILE A 138	-62.770	45.131	12.238	1.00	30.20	C
ATOM	1828	CG1	ILE A 138	-62.914	46.576	12.689	1.00	29.92	C
ATOM	1831	CD1	ILE A 138	-61.681	47.414	12.417	1.00	30.50	C
ATOM	1835	CG2	ILE A 138	-61.706	44.470	13.085	1.00	29.73	C
ATOM	1839	C	ILE A 138	-62.140	43.631	10.320	1.00	30.34	C
ATOM	1840	O	ILE A 138	-60.972	43.203	10.212	1.00	30.65	O
ATOM	1842	N	LEU A 139	-63.227	42.884	10.098	1.00	29.69	N
ATOM	1843	CA	LEU A 139	-63.134	41.523	9.562	1.00	29.23	C
ATOM	1845	CB	LEU A 139	-64.517	40.917	9.332	1.00	28.92	C
ATOM	1848	CG	LEU A 139	-65.066	40.185	10.545	1.00	29.06	C
ATOM	1850	CD1	LEU A 139	-66.541	39.875	10.352	1.00	28.20	C
ATOM	1854	CD2	LEU A 139	-64.250	38.922	10.813	1.00	28.59	C
ATOM	1858	C	LEU A 139	-62.382	41.546	8.247	1.00	29.00	C
ATOM	1859	O	LEU A 139	-61.444	40.773	8.035	1.00	29.50	O

ATOM	1861	N	SER A 140	-62.790	42.458	7.372	1.00	28.14	N
ATOM	1862	CA	SER A 140	-62.222	42.526	6.054	1.00	27.40	C
ATOM	1864	CB	SER A 140	-62.912	43.619	5.260	1.00	27.35	C
ATOM	1867	OG	SER A 140	-63.310	43.103	4.015	1.00	28.43	O
ATOM	1869	C	SER A 140	-60.720	42.777	6.123	1.00	26.77	C
ATOM	1870	O	SER A 140	-59.975	42.269	5.298	1.00	27.29	O
ATOM	1872	N	LEU A 141	-60.285	43.563	7.105	1.00	25.89	N
ATOM	1873	CA	LEU A 141	-58.872	43.897	7.283	1.00	25.20	C
ATOM	1875	CB	LEU A 141	-58.717	45.164	8.155	1.00	25.00	C
ATOM	1878	CG	LEU A 141	-57.298	45.670	8.489	1.00	24.35	C
ATOM	1880	CD1	LEU A 141	-56.497	45.974	7.232	1.00	22.75	C
ATOM	1884	CD2	LEU A 141	-57.329	46.893	9.385	1.00	22.81	C
ATOM	1888	C	LEU A 141	-58.140	42.721	7.920	1.00	24.92	C
ATOM	1889	O	LEU A 141	-57.034	42.380	7.521	1.00	24.77	O
ATOM	1891	N	TYR A 142	-58.751	42.101	8.921	1.00	24.51	N
ATOM	1892	CA	TYR A 142	-58.163	40.913	9.514	1.00	24.36	C
ATOM	1894	CB	TYR A 142	-59.120	40.325	10.538	1.00	24.18	C
ATOM	1897	CG	TYR A 142	-58.774	38.941	11.046	1.00	23.72	C
ATOM	1898	CD1	TYR A 142	-57.878	38.763	12.091	1.00	22.97	C
ATOM	1900	CE1	TYR A 142	-57.593	37.500	12.582	1.00	24.25	C
ATOM	1902	CZ	TYR A 142	-58.229	36.386	12.030	1.00	25.75	C
ATOM	1903	OH	TYR A 142	-57.967	35.100	12.500	1.00	26.97	O
ATOM	1905	CE2	TYR A 142	-59.120	36.551	10.988	1.00	25.02	C
ATOM	1907	CD2	TYR A 142	-59.390	37.820	10.511	1.00	24.01	C
ATOM	1909	C	TYR A 142	-57.877	39.896	8.423	1.00	24.39	C
ATOM	1910	O	TYR A 142	-56.822	39.276	8.380	1.00	24.19	O
ATOM	1912	N	GLU A 143	-58.841	39.760	7.527	1.00	24.59	N
ATOM	1913	CA	GLU A 143	-58.820	38.728	6.521	1.00	24.70	C
ATOM	1915	CB	GLU A 143	-60.194	38.636	5.853	1.00	24.69	C
ATOM	1918	CG	GLU A 143	-60.651	37.219	5.547	1.00	27.21	C
ATOM	1921	CD	GLU A 143	-61.127	36.426	6.774	1.00	30.20	C
ATOM	1922	OE1	GLU A 143	-60.261	35.776	7.393	1.00	33.40	O
ATOM	1923	OE2	GLU A 143	-62.351	36.421	7.097	1.00	30.50	O
ATOM	1924	C	GLU A 143	-57.701	39.010	5.518	1.00	24.08	C
ATOM	1925	O	GLU A 143	-56.984	38.090	5.126	1.00	24.53	O
ATOM	1927	N	ALA A 144	-57.537	40.281	5.147	1.00	23.30	N
ATOM	1928	CA	ALA A 144	-56.515	40.715	4.173	1.00	22.79	C
ATOM	1930	CB	ALA A 144	-56.787	42.154	3.740	1.00	22.49	C
ATOM	1934	C	ALA A 144	-55.073	40.598	4.701	1.00	22.42	C
ATOM	1935	O	ALA A 144	-54.128	40.380	3.935	1.00	22.36	O
ATOM	1937	N	SER A 145	-54.909	40.713	6.012	1.00	21.80	N
ATOM	1938	CA	SER A 145	-53.593	40.822	6.600	1.00	21.36	C
ATOM	1940	CB	SER A 145	-53.726	41.170	8.066	1.00	21.21	C
ATOM	1943	OG	SER A 145	-54.337	40.100	8.753	1.00	21.59	O
ATOM	1945	C	SER A 145	-52.797	39.537	6.476	1.00	21.30	C
ATOM	1946	O	SER A 145	-51.571	39.533	6.701	1.00	21.50	O
ATOM	1948	N	PHE A 146	-53.488	38.439	6.154	1.00	21.00	N
ATOM	1949	CA	PHE A 146	-52.827	37.127	6.030	1.00	20.41	C
ATOM	1951	CB	PHE A 146	-53.796	35.979	6.338	1.00	20.20	C
ATOM	1954	CG	PHE A 146	-54.130	35.873	7.798	1.00	19.91	C
ATOM	1955	CD1	PHE A 146	-53.382	35.087	8.637	1.00	20.48	C
ATOM	1957	CE1	PHE A 146	-53.679	35.012	9.992	1.00	20.59	C
ATOM	1959	CZ	PHE A 146	-54.731	35.740	10.508	1.00	19.15	C
ATOM	1961	CE2	PHE A 146	-55.464	36.536	9.691	1.00	19.10	C
ATOM	1963	CD2	PHE A 146	-55.155	36.615	8.343	1.00	20.16	C
ATOM	1965	C	PHE A 146	-52.174	36.964	4.680	1.00	19.76	C
ATOM	1966	O	PHE A 146	-51.305	36.116	4.523	1.00	19.46	O
ATOM	1968	N	LEU A 147	-52.550	37.819	3.729	1.00	19.34	N
ATOM	1969	CA	LEU A 147	-51.933	37.804	2.402	1.00	19.13	C
ATOM	1971	CB	LEU A 147	-52.905	38.367	1.347	1.00	18.67	C

ATOM	1974	CG	LEU	A	147	-53.964	37.334	.919	1.00	18.11	C
ATOM	1976	CD1	LEU	A	147	-55.090	37.248	1.961	1.00	15.10	C
ATOM	1980	CD2	LEU	A	147	-54.494	37.609	-.499	1.00	16.50	C
ATOM	1984	C	LEU	A	147	-50.582	38.525	2.369	1.00	19.10	C
ATOM	1985	O	LEU	A	147	-49.992	38.715	1.311	1.00	19.44	O
ATOM	1987	N	ALA	A	148	-50.075	38.882	3.542	1.00	19.13	N
ATOM	1988	CA	ALA	A	148	-48.967	39.806	3.664	1.00	18.84	C
ATOM	1990	CB	ALA	A	148	-48.762	40.167	5.117	1.00	18.78	C
ATOM	1994	C	ALA	A	148	-47.723	39.184	3.113	1.00	18.85	C
ATOM	1995	O	ALA	A	148	-47.548	37.975	3.206	1.00	18.80	O
ATOM	1997	N	LEU	A	149	-46.869	40.013	2.523	1.00	19.08	N
ATOM	1998	CA	LEU	A	149	-45.507	39.601	2.215	1.00	19.32	C
ATOM	2000	CB	LEU	A	149	-45.056	40.145	.863	1.00	19.14	C
ATOM	2003	CG	LEU	A	149	-45.819	39.650	-.364	1.00	19.16	C
ATOM	2005	CD1	LEU	A	149	-45.036	39.929	-1.645	1.00	18.43	C
ATOM	2009	CD2	LEU	A	149	-46.089	38.175	-.259	1.00	20.09	C
ATOM	2013	C	LEU	A	149	-44.527	40.029	3.304	1.00	19.60	C
ATOM	2014	O	LEU	A	149	-44.797	40.916	4.113	1.00	19.20	O
ATOM	2016	N	GLU	A	150	-43.384	39.361	3.309	1.00	20.42	N
ATOM	2017	CA	GLU	A	150	-42.286	39.674	4.215	1.00	21.03	C
ATOM	2019	CB	GLU	A	150	-41.098	38.741	3.948	1.00	21.03	C
ATOM	2022	CG	GLU	A	150	-40.491	38.176	5.205	1.00	22.23	C
ATOM	2025	CD	GLU	A	150	-39.433	37.117	4.941	1.00	24.34	C
ATOM	2026	OE1	GLU	A	150	-39.189	36.786	3.753	1.00	24.53	O
ATOM	2027	OE2	GLU	A	150	-38.843	36.625	5.940	1.00	25.82	O
ATOM	2028	C	GLU	A	150	-41.887	41.130	4.011	1.00	21.29	C
ATOM	2029	O	GLU	A	150	-41.619	41.544	2.885	1.00	21.13	O
ATOM	2031	N	GLY	A	151	-41.893	41.903	5.094	1.00	21.89	N
ATOM	2032	CA	GLY	A	151	-41.636	43.341	5.030	1.00	22.25	C
ATOM	2035	C	GLY	A	151	-42.871	44.221	5.152	1.00	22.67	C
ATOM	2036	O	GLY	A	151	-42.757	45.397	5.462	1.00	22.84	O
ATOM	2038	N	GLU	A	152	-44.055	43.669	4.918	1.00	23.19	N
ATOM	2039	CA	GLU	A	152	-45.263	44.474	4.937	1.00	23.76	C
ATOM	2041	CB	GLU	A	152	-46.350	43.813	4.100	1.00	23.79	C
ATOM	2044	CG	GLU	A	152	-46.037	43.890	2.615	1.00	23.98	C
ATOM	2047	CD	GLU	A	152	-47.147	43.368	1.749	1.00	23.72	C
ATOM	2048	OE1	GLU	A	152	-47.226	43.771	.574	1.00	22.11	O
ATOM	2049	OE2	GLU	A	152	-47.939	42.542	2.241	1.00	25.41	O
ATOM	2050	C	GLU	A	152	-45.722	44.719	6.365	1.00	24.47	C
ATOM	2051	O	GLU	A	152	-46.586	44.006	6.906	1.00	24.70	O
ATOM	2053	N	ASN	A	153	-45.134	45.744	6.969	1.00	25.00	N
ATOM	2054	CA	ASN	A	153	-45.287	45.978	8.404	1.00	25.59	C
ATOM	2056	CB	ASN	A	153	-44.221	46.957	8.902	1.00	25.55	C
ATOM	2059	CG	ASN	A	153	-44.386	48.331	8.308	1.00	25.57	C
ATOM	2060	OD1	ASN	A	153	-44.114	48.557	7.128	1.00	24.52	O
ATOM	2061	ND2	ASN	A	153	-44.867	49.256	9.117	1.00	27.05	N
ATOM	2064	C	ASN	A	153	-46.679	46.486	8.772	1.00	25.92	C
ATOM	2065	O	ASN	A	153	-47.172	46.238	9.872	1.00	26.02	O
ATOM	2067	N	ILE	A	154	-47.310	47.197	7.847	1.00	26.35	N
ATOM	2068	CA	ILE	A	154	-48.638	47.751	8.082	1.00	26.57	C
ATOM	2070	CB	ILE	A	154	-49.013	48.779	7.007	1.00	26.45	C
ATOM	2072	CG1	ILE	A	154	-48.159	50.023	7.164	1.00	26.36	C
ATOM	2075	CD1	ILE	A	154	-48.060	50.790	5.883	1.00	27.92	C
ATOM	2079	CG2	ILE	A	154	-50.467	49.152	7.089	1.00	26.25	C
ATOM	2083	C	ILE	A	154	-49.674	46.643	8.111	1.00	26.96	C
ATOM	2084	O	ILE	A	154	-50.650	46.741	8.851	1.00	27.12	O
ATOM	2086	N	LEU	A	155	-49.467	45.588	7.320	1.00	27.35	N
ATOM	2087	CA	LEU	A	155	-50.422	44.475	7.304	1.00	27.76	C
ATOM	2089	CB	LEU	A	155	-50.238	43.594	6.067	1.00	27.73	C
ATOM	2092	CG	LEU	A	155	-51.090	43.990	4.875	1.00	27.31	C

ATOM	2094	CD1 LEU A 155	-50.980	45.456	4.683	1.00	28.13	C
ATOM	2098	CD2 LEU A 155	-50.633	43.262	3.623	1.00	27.79	C
ATOM	2102	C LEU A 155	-50.340	43.642	8.582	1.00	28.22	C
ATOM	2103	O LEU A 155	-51.361	43.237	9.112	1.00	27.95	O
ATOM	2105	N ASP A 156	-49.128	43.395	9.068	1.00	28.95	N
ATOM	2106	CA ASP A 156	-48.944	42.723	10.348	1.00	29.71	C
ATOM	2108	CB ASP A 156	-47.450	42.461	10.639	1.00	30.23	C
ATOM	2111	CG ASP A 156	-46.938	41.125	10.027	1.00	32.82	C
ATOM	2112	OD1 ASP A 156	-47.765	40.240	9.658	1.00	36.27	O
ATOM	2113	OD2 ASP A 156	-45.700	40.961	9.923	1.00	34.48	O
ATOM	2114	C ASP A 156	-49.570	43.539	11.469	1.00	29.66	C
ATOM	2115	O ASP A 156	-50.109	42.969	12.425	1.00	29.81	O
ATOM	2117	N GLU A 157	-49.514	44.866	11.345	1.00	29.69	N
ATOM	2118	CA GLU A 157	-50.129	45.768	12.329	1.00	29.57	C
ATOM	2120	CB GLU A 157	-49.581	47.186	12.178	1.00	29.71	C
ATOM	2123	CG GLU A 157	-48.164	47.326	12.712	1.00	30.53	C
ATOM	2126	CD GLU A 157	-47.455	48.544	12.169	1.00	31.96	C
ATOM	2127	OE1 GLU A 157	-48.167	49.491	11.741	1.00	33.56	O
ATOM	2128	OE2 GLU A 157	-46.195	48.546	12.166	1.00	31.54	O
ATOM	2129	C GLU A 157	-51.649	45.771	12.230	1.00	29.16	C
ATOM	2130	O GLU A 157	-52.330	45.819	13.241	1.00	29.00	O
ATOM	2132	N ALA A 158	-52.172	45.716	11.010	1.00	28.86	N
ATOM	2133	CA ALA A 158	-53.606	45.570	10.792	1.00	28.79	C
ATOM	2135	CB ALA A 158	-53.893	45.456	9.301	1.00	28.44	C
ATOM	2139	C ALA A 158	-54.171	44.352	11.536	1.00	28.93	C
ATOM	2140	O ALA A 158	-55.265	44.402	12.098	1.00	28.77	O
ATOM	2142	N LYS A 159	-53.416	43.260	11.528	1.00	29.19	N
ATOM	2143	CA LYS A 159	-53.849	42.022	12.143	1.00	29.57	C
ATOM	2145	CB LYS A 159	-52.929	40.880	11.710	1.00	29.62	C
ATOM	2148	CG LYS A 159	-53.297	39.513	12.258	1.00	30.27	C
ATOM	2151	CD LYS A 159	-52.765	38.365	11.387	1.00	31.78	C
ATOM	2154	CE LYS A 159	-51.259	38.132	11.525	1.00	33.15	C
ATOM	2157	NZ LYS A 159	-50.911	36.696	11.248	1.00	33.88	N
ATOM	2161	C LYS A 159	-53.866	42.160	13.670	1.00	29.99	C
ATOM	2162	O LYS A 159	-54.767	41.642	14.332	1.00	30.29	O
ATOM	2164	N VAL A 160	-52.886	42.862	14.236	1.00	30.07	N
ATOM	2165	CA VAL A 160	-52.840	43.047	15.696	1.00	30.03	C
ATOM	2167	CB VAL A 160	-51.520	43.691	16.162	1.00	29.88	C
ATOM	2169	CG1 VAL A 160	-51.515	43.863	17.671	1.00	28.94	C
ATOM	2173	CG2 VAL A 160	-50.341	42.848	15.707	1.00	30.23	C
ATOM	2177	C VAL A 160	-53.996	43.925	16.149	1.00	30.05	C
ATOM	2178	O VAL A 160	-54.544	43.752	17.234	1.00	30.31	O
ATOM	2180	N PHE A 161	-54.355	44.862	15.290	1.00	30.08	N
ATOM	2181	CA PHE A 161	-55.407	45.808	15.558	1.00	30.09	C
ATOM	2183	CB PHE A 161	-55.255	46.993	14.598	1.00	29.92	C
ATOM	2186	CG PHE A 161	-56.370	47.974	14.654	1.00	28.52	C
ATOM	2187	CD1 PHE A 161	-56.347	49.005	15.554	1.00	27.92	C
ATOM	2189	CE1 PHE A 161	-57.361	49.903	15.600	1.00	28.06	C
ATOM	2191	CZ PHE A 161	-58.411	49.783	14.733	1.00	28.76	C
ATOM	2193	CE2 PHE A 161	-58.441	48.757	13.826	1.00	28.31	C
ATOM	2195	CD2 PHE A 161	-57.429	47.869	13.789	1.00	28.16	C
ATOM	2197	C PHE A 161	-56.746	45.123	15.388	1.00	30.55	C
ATOM	2198	O PHE A 161	-57.598	45.203	16.258	1.00	30.53	O
ATOM	2200	N ALA A 162	-56.930	44.438	14.271	1.00	31.34	N
ATOM	2201	CA ALA A 162	-58.193	43.770	14.023	1.00	32.26	C
ATOM	2203	CB ALA A 162	-58.217	43.139	12.643	1.00	32.31	C
ATOM	2207	C ALA A 162	-58.500	42.734	15.118	1.00	33.09	C
ATOM	2208	O ALA A 162	-59.557	42.814	15.735	1.00	33.83	O
ATOM	2210	N ILE A 163	-57.587	41.799	15.393	1.00	33.59	N
ATOM	2211	CA ILE A 163	-57.826	40.799	16.438	1.00	34.02	C

ATOM	2213	CB	ILE A 163	-56.597	39.907	16.745	1.00	33.98	C
ATOM	2215	CG1	ILE A 163	-56.235	39.008	15.566	1.00	34.37	C
ATOM	2218	CD1	ILE A 163	-54.878	38.317	15.722	1.00	34.38	C
ATOM	2222	CG2	ILE A 163	-56.883	38.995	17.919	1.00	33.38	C
ATOM	2226	C	ILE A 163	-58.225	41.472	17.744	1.00	34.75	C
ATOM	2227	O	ILE A 163	-59.189	41.064	18.373	1.00	35.03	O
ATOM	2229	N	SER A 164	-57.501	42.504	18.156	1.00	35.76	N
ATOM	2230	CA	SER A 164	-57.763	43.117	19.463	1.00	36.59	C
ATOM	2232	CB	SER A 164	-56.828	44.293	19.735	1.00	36.59	C
ATOM	2235	OG	SER A 164	-57.350	45.485	19.164	1.00	36.46	O
ATOM	2237	C	SER A 164	-59.196	43.611	19.556	1.00	37.42	C
ATOM	2238	O	SER A 164	-59.860	43.420	20.571	1.00	37.38	O
ATOM	2240	N	HIS A 165	-59.661	44.251	18.486	1.00	38.51	N
ATOM	2241	CA	HIS A 165	-60.990	44.851	18.467	1.00	39.40	C
ATOM	2243	CB	HIS A 165	-61.008	46.097	17.574	1.00	39.76	C
ATOM	2246	CG	HIS A 165	-60.467	47.323	18.251	1.00	42.06	C
ATOM	2247	ND1	HIS A 165	-59.205	47.821	18.000	1.00	44.21	N
ATOM	2249	CE1	HIS A 165	-58.998	48.893	18.748	1.00	45.04	C
ATOM	2251	NE2	HIS A 165	-60.080	49.106	19.478	1.00	44.99	N
ATOM	2253	CD2	HIS A 165	-61.011	48.135	19.192	1.00	43.76	C
ATOM	2255	C	HIS A 165	-62.078	43.872	18.060	1.00	39.38	C
ATOM	2256	O	HIS A 165	-63.248	44.244	18.067	1.00	39.43	O
ATOM	2258	N	LEU A 166	-61.683	42.632	17.745	1.00	39.54	N
ATOM	2259	CA	LEU A 166	-62.589	41.550	17.326	1.00	39.55	C
ATOM	2261	CB	LEU A 166	-62.001	40.766	16.154	1.00	39.28	C
ATOM	2264	CG	LEU A 166	-62.156	41.339	14.752	1.00	38.23	C
ATOM	2266	CD1	LEU A 166	-61.361	40.501	13.805	1.00	37.53	C
ATOM	2270	CD2	LEU A 166	-63.594	41.367	14.332	1.00	37.21	C
ATOM	2274	C	LEU A 166	-62.842	40.555	18.438	1.00	40.08	C
ATOM	2275	O	LEU A 166	-63.991	40.220	18.733	1.00	40.42	O
ATOM	2277	N	LYS A 167	-61.766	40.074	19.052	1.00	40.58	N
ATOM	2278	CA	LYS A 167	-61.875	39.090	20.133	1.00	41.09	C
ATOM	2280	CB	LYS A 167	-60.491	38.632	20.595	1.00	41.13	C
ATOM	2283	CG	LYS A 167	-59.897	39.399	21.768	1.00	42.16	C
ATOM	2286	CD	LYS A 167	-58.421	39.035	22.013	1.00	44.45	C
ATOM	2289	CE	LYS A 167	-58.070	37.561	21.654	1.00	45.56	C
ATOM	2292	NZ	LYS A 167	-56.610	37.237	21.798	1.00	46.20	N
ATOM	2296	C	LYS A 167	-62.678	39.584	21.337	1.00	41.23	C
ATOM	2297	O	LYS A 167	-62.994	38.793	22.215	1.00	41.47	O
ATOM	2299	N	GLU A 168	-62.995	40.884	21.352	1.00	41.49	N
ATOM	2300	CA	GLU A 168	-63.704	41.583	22.437	1.00	41.42	C
ATOM	2302	CB	GLU A 168	-65.201	41.730	22.100	1.00	41.27	C
ATOM	2305	CG	GLU A 168	-65.995	40.424	22.016	1.00	41.01	C
ATOM	2308	CD	GLU A 168	-66.928	40.328	20.787	1.00	40.84	C
ATOM	2309	OE1	GLU A 168	-66.701	41.040	19.770	1.00	39.82	O
ATOM	2310	OE2	GLU A 168	-67.884	39.509	20.845	1.00	39.04	O
ATOM	2311	C	GLU A 168	-63.484	41.004	23.843	1.00	41.66	C
ATOM	2312	O	GLU A 168	-62.435	41.231	24.469	1.00	41.46	O
ATOM	2314	N	GLU A 172	-71.668	38.961	21.192	1.00	55.65	N
ATOM	2315	CA	GLU A 172	-73.052	39.153	20.721	1.00	56.00	C
ATOM	2317	CB	GLU A 172	-73.876	37.880	20.937	1.00	56.27	C
ATOM	2320	CG	GLU A 172	-73.198	36.620	20.390	1.00	57.12	C
ATOM	2323	CD	GLU A 172	-74.056	35.363	20.516	1.00	58.18	C
ATOM	2324	OE1	GLU A 172	-75.300	35.480	20.622	1.00	58.84	O
ATOM	2325	OE2	GLU A 172	-73.480	34.249	20.498	1.00	58.98	O
ATOM	2326	C	GLU A 172	-73.766	40.386	21.334	1.00	55.70	C
ATOM	2327	O	GLU A 172	-74.956	40.347	21.688	1.00	54.95	O
ATOM	2329	N	LYS A 173	-72.987	41.458	21.486	1.00	55.58	N
ATOM	2330	CA	LYS A 173	-73.496	42.831	21.466	1.00	55.43	C
ATOM	2332	CB	LYS A 173	-72.649	43.741	22.357	1.00	55.51	C

ATOM	2335	CG	LYS A 173	-72.527	43.223	23.793	1.00	56.07	C
ATOM	2338	CD	LYS A 173	-72.803	44.291	24.868	1.00	55.87	C
ATOM	2341	CE	LYS A 173	-73.111	43.642	26.221	1.00	55.20	C
ATOM	2344	NZ	LYS A 173	-72.893	44.576	27.348	1.00	54.63	N
ATOM	2348	C	LYS A 173	-73.449	43.297	20.009	1.00	55.01	C
ATOM	2349	O	LYS A 173	-74.091	44.274	19.625	1.00	54.97	O
ATOM	2351	N	ILE A 174	-72.647	42.581	19.223	1.00	54.57	N
ATOM	2352	CA	ILE A 174	-72.713	42.550	17.767	1.00	54.22	C
ATOM	2354	CB	ILE A 174	-71.429	41.845	17.227	1.00	54.21	C
ATOM	2356	CG1	ILE A 174	-70.196	42.724	17.474	1.00	54.73	C
ATOM	2359	CD1	ILE A 174	-68.858	41.949	17.553	1.00	55.78	C
ATOM	2363	CG2	ILE A 174	-71.532	41.504	15.759	1.00	54.06	C
ATOM	2367	C	ILE A 174	-73.957	41.744	17.384	1.00	53.88	C
ATOM	2368	O	ILE A 174	-74.587	41.134	18.250	1.00	53.83	O
ATOM	2370	N	GLY A 175	-74.332	41.747	16.107	1.00	53.59	N
ATOM	2371	CA	GLY A 175	-75.351	40.817	15.613	1.00	53.51	C
ATOM	2374	C	GLY A 175	-75.130	39.395	16.126	1.00	53.43	C
ATOM	2375	O	GLY A 175	-74.119	39.105	16.760	1.00	53.80	O
ATOM	2377	N	LYS A 176	-76.070	38.495	15.861	1.00	53.10	N
ATOM	2378	CA	LYS A 176	-75.926	37.098	16.286	1.00	52.68	C
ATOM	2380	CB	LYS A 176	-77.291	36.522	16.677	1.00	52.99	C
ATOM	2383	CG	LYS A 176	-78.003	37.387	17.746	1.00	54.12	C
ATOM	2386	CD	LYS A 176	-79.201	36.697	18.436	1.00	54.57	C
ATOM	2389	CE	LYS A 176	-79.503	37.369	19.792	1.00	54.79	C
ATOM	2392	NZ	LYS A 176	-80.458	36.598	20.653	1.00	55.20	N
ATOM	2396	C	LYS A 176	-75.246	36.267	15.196	1.00	51.83	C
ATOM	2397	O	LYS A 176	-74.393	35.433	15.486	1.00	51.69	O
ATOM	2399	N	GLU A 177	-75.609	36.520	13.940	1.00	50.90	N
ATOM	2400	CA	GLU A 177	-74.975	35.854	12.800	1.00	50.19	C
ATOM	2402	CB	GLU A 177	-75.792	36.062	11.507	1.00	50.41	C
ATOM	2405	CG	GLU A 177	-75.506	37.368	10.735	1.00	51.03	C
ATOM	2408	CD	GLU A 177	-76.403	37.572	9.507	1.00	52.17	C
ATOM	2409	OE1	GLU A 177	-75.955	38.276	8.569	1.00	52.63	O
ATOM	2410	OE2	GLU A 177	-77.546	37.048	9.479	1.00	52.03	O
ATOM	2411	C	GLU A 177	-73.549	36.369	12.607	1.00	49.10	C
ATOM	2412	O	GLU A 177	-72.626	35.597	12.321	1.00	49.53	O
ATOM	2414	N	LEU A 178	-73.388	37.681	12.770	1.00	47.39	N
ATOM	2415	CA	LEU A 178	-72.110	38.371	12.575	1.00	45.63	C
ATOM	2417	CB	LEU A 178	-72.352	39.877	12.615	1.00	45.59	C
ATOM	2420	CG	LEU A 178	-71.280	40.809	12.085	1.00	45.45	C
ATOM	2422	CD1	LEU A 178	-70.969	40.515	10.620	1.00	45.45	C
ATOM	2426	CD2	LEU A 178	-71.765	42.233	12.279	1.00	44.63	C
ATOM	2430	C	LEU A 178	-71.095	37.979	13.634	1.00	44.06	C
ATOM	2431	O	LEU A 178	-69.902	38.012	13.387	1.00	43.70	O
ATOM	2433	N	ALA A 179	-71.586	37.611	14.810	1.00	42.54	N
ATOM	2434	CA	ALA A 179	-70.753	37.041	15.847	1.00	41.66	C
ATOM	2436	CB	ALA A 179	-71.545	36.884	17.122	1.00	41.69	C
ATOM	2440	C	ALA A 179	-70.193	35.695	15.416	1.00	40.75	C
ATOM	2441	O	ALA A 179	-69.078	35.351	15.769	1.00	40.61	O
ATOM	2443	N	GLU A 180	-70.974	34.931	14.664	1.00	39.88	N
ATOM	2444	CA	GLU A 180	-70.518	33.639	14.143	1.00	39.40	C
ATOM	2446	CB	GLU A 180	-71.710	32.781	13.718	1.00	39.74	C
ATOM	2449	CG	GLU A 180	-72.650	32.415	14.865	1.00	40.88	C
ATOM	2452	CD	GLU A 180	-73.913	31.710	14.398	1.00	42.02	C
ATOM	2453	OE1	GLU A 180	-74.380	32.022	13.271	1.00	43.05	O
ATOM	2454	OE2	GLU A 180	-74.437	30.857	15.166	1.00	41.93	O
ATOM	2455	C	GLU A 180	-69.565	33.805	12.960	1.00	38.29	C
ATOM	2456	O	GLU A 180	-68.739	32.940	12.699	1.00	38.23	O
ATOM	2458	N	GLN A 181	-69.704	34.911	12.240	1.00	36.92	N
ATOM	2459	CA	GLN A 181	-68.748	35.299	11.223	1.00	35.74	C

ATOM	2461	CB	GLN A 181	-69.272	36.516	10.463	1.00	35.94	C
ATOM	2464	CG	GLN A 181	-68.931	36.537	8.975	1.00	36.93	C
ATOM	2467	CD	GLN A 181	-69.693	35.492	8.176	1.00	37.72	C
ATOM	2468	OE1	GLN A 181	-70.769	35.067	8.567	1.00	38.78	O
ATOM	2469	NE2	GLN A 181	-69.139	35.083	7.050	1.00	38.11	N
ATOM	2472	C	GLN A 181	-67.406	35.632	11.881	1.00	34.64	C
ATOM	2473	O	GLN A 181	-66.348	35.411	11.294	1.00	34.73	O
ATOM	2475	N	VAL A 182	-67.448	36.158	13.102	1.00	33.29	N
ATOM	2476	CA	VAL A 182	-66.224	36.488	13.842	1.00	32.28	C
ATOM	2478	CB	VAL A 182	-66.433	37.653	14.855	1.00	32.10	C
ATOM	2480	CG1	VAL A 182	-66.328	38.991	14.166	1.00	31.51	C
ATOM	2484	CG2	VAL A 182	-65.420	37.594	15.957	1.00	31.69	C
ATOM	2488	C	VAL A 182	-65.613	35.298	14.575	1.00	31.49	C
ATOM	2489	O	VAL A 182	-64.396	35.230	14.668	1.00	31.43	O
ATOM	2491	N	SER A 183	-66.422	34.378	15.106	1.00	30.56	N
ATOM	2492	CA	SER A 183	-65.860	33.174	15.744	1.00	30.16	C
ATOM	2494	CB	SER A 183	-66.926	32.253	16.327	1.00	30.29	C
ATOM	2497	OG	SER A 183	-67.796	32.934	17.214	1.00	32.35	O
ATOM	2499	C	SER A 183	-65.111	32.403	14.692	1.00	29.30	C
ATOM	2500	O	SER A 183	-63.974	31.981	14.903	1.00	29.32	O
ATOM	2502	N	HIS A 184	-65.768	32.252	13.546	1.00	28.24	N
ATOM	2503	CA	HIS A 184	-65.215	31.578	12.379	1.00	27.36	C
ATOM	2505	CB	HIS A 184	-66.245	31.653	11.254	1.00	27.41	C
ATOM	2508	CG	HIS A 184	-65.930	30.807	10.064	1.00	26.90	C
ATOM	2509	ND1	HIS A 184	-65.484	29.511	10.166	1.00	27.76	N
ATOM	2511	CE1	HIS A 184	-65.321	29.009	8.955	1.00	27.47	C
ATOM	2513	NE2	HIS A 184	-65.661	29.929	8.075	1.00	26.52	N
ATOM	2515	CD2	HIS A 184	-66.059	31.058	8.742	1.00	26.61	C
ATOM	2517	C	HIS A 184	-63.898	32.200	11.940	1.00	26.69	C
ATOM	2518	O	HIS A 184	-62.881	31.525	11.862	1.00	26.44	O
ATOM	2520	N	ALA A 185	-63.909	33.498	11.675	1.00	26.11	N
ATOM	2521	CA	ALA A 185	-62.711	34.176	11.205	1.00	25.63	C
ATOM	2523	CB	ALA A 185	-62.970	35.627	11.109	1.00	25.74	C
ATOM	2527	C	ALA A 185	-61.552	33.917	12.153	1.00	25.33	C
ATOM	2528	O	ALA A 185	-60.451	33.572	11.725	1.00	25.79	O
ATOM	2530	N	LEU A 186	-61.814	34.066	13.450	1.00	24.77	N
ATOM	2531	CA	LEU A 186	-60.776	33.940	14.470	1.00	24.05	C
ATOM	2533	CB	LEU A 186	-61.256	34.498	15.814	1.00	23.63	C
ATOM	2536	CG	LEU A 186	-61.516	36.012	15.903	1.00	22.51	C
ATOM	2538	CD1	LEU A 186	-62.057	36.399	17.266	1.00	20.79	C
ATOM	2542	CD2	LEU A 186	-60.273	36.793	15.613	1.00	21.13	C
ATOM	2546	C	LEU A 186	-60.316	32.497	14.613	1.00	24.02	C
ATOM	2547	O	LEU A 186	-59.169	32.250	14.943	1.00	24.27	O
ATOM	2549	N	GLU A 187	-61.201	31.543	14.347	1.00	23.88	N
ATOM	2550	CA	GLU A 187	-60.836	30.117	14.348	1.00	23.53	C
ATOM	2552	CB	GLU A 187	-62.029	29.270	13.919	1.00	23.66	C
ATOM	2555	CG	GLU A 187	-61.945	27.828	14.316	1.00	25.26	C
ATOM	2558	CD	GLU A 187	-63.049	26.971	13.693	1.00	28.04	C
ATOM	2559	OE1	GLU A 187	-63.672	27.393	12.681	1.00	28.67	O
ATOM	2560	OE2	GLU A 187	-63.288	25.858	14.225	1.00	29.72	O
ATOM	2561	C	GLU A 187	-59.686	29.838	13.395	1.00	22.92	C
ATOM	2562	O	GLU A 187	-58.765	29.091	13.729	1.00	22.53	O
ATOM	2564	N	LEU A 188	-59.763	30.446	12.209	1.00	22.46	N
ATOM	2565	CA	LEU A 188	-58.799	30.246	11.121	1.00	22.07	C
ATOM	2567	CB	LEU A 188	-58.802	28.794	10.635	1.00	22.11	C
ATOM	2570	CG	LEU A 188	-57.876	28.362	9.493	1.00	22.17	C
ATOM	2572	CD1	LEU A 188	-56.456	28.051	9.961	1.00	21.62	C
ATOM	2576	CD2	LEU A 188	-58.464	27.134	8.822	1.00	22.40	C
ATOM	2580	C	LEU A 188	-59.237	31.137	9.985	1.00	21.69	C
ATOM	2581	O	LEU A 188	-60.380	31.075	9.563	1.00	21.53	O

ATOM	2583	N	PRO A 189	-58.333	31.973	9.474	1.00	21.46	N
ATOM	2584	CA	PRO A 189	-58.745	32.935	8.472	1.00	21.18	C
ATOM	2586	CB	PRO A 189	-57.592	33.936	8.486	1.00	21.03	C
ATOM	2589	CG	PRO A 189	-56.423	33.118	8.757	1.00	21.21	C
ATOM	2592	CD	PRO A 189	-56.873	31.992	9.668	1.00	21.61	C
ATOM	2595	C	PRO A 189	-58.946	32.314	7.076	1.00	20.96	C
ATOM	2596	O	PRO A 189	-58.418	31.247	6.793	1.00	21.02	O
ATOM	2597	N	LEU A 190	-59.714	32.994	6.226	1.00	20.89	N
ATOM	2598	CA	LEU A 190	-60.018	32.533	4.867	1.00	20.74	C
ATOM	2600	CB	LEU A 190	-60.597	33.690	4.053	1.00	20.76	C
ATOM	2603	CG	LEU A 190	-62.090	33.924	4.259	1.00	21.35	C
ATOM	2605	CD1	LEU A 190	-62.537	35.196	3.578	1.00	20.54	C
ATOM	2609	CD2	LEU A 190	-62.889	32.712	3.720	1.00	23.39	C
ATOM	2613	C	LEU A 190	-58.814	31.971	4.112	1.00	20.35	C
ATOM	2614	O	LEU A 190	-58.859	30.871	3.546	1.00	20.08	O
ATOM	2616	N	HIS A 191	-57.734	32.746	4.132	1.00	19.81	N
ATOM	2617	CA	HIS A 191	-56.537	32.448	3.364	1.00	18.89	C
ATOM	2619	CB	HIS A 191	-55.548	33.592	3.536	1.00	18.92	C
ATOM	2622	CG	HIS A 191	-54.372	33.520	2.623	1.00	18.51	C
ATOM	2623	ND1	HIS A 191	-54.482	33.647	1.258	1.00	18.08	N
ATOM	2625	CE1	HIS A 191	-53.281	33.544	.718	1.00	18.96	C
ATOM	2627	NE2	HIS A 191	-52.397	33.371	1.684	1.00	17.01	N
ATOM	2629	CD2	HIS A 191	-53.054	33.355	2.884	1.00	17.94	C
ATOM	2631	C	HIS A 191	-55.889	31.142	3.759	1.00	18.32	C
ATOM	2632	O	HIS A 191	-55.217	30.531	2.942	1.00	17.98	O
ATOM	2634	N	ARG A 192	-56.093	30.719	5.006	1.00	18.06	N
ATOM	2635	CA	ARG A 192	-55.502	29.480	5.523	1.00	17.99	C
ATOM	2637	CB	ARG A 192	-54.940	29.722	6.909	1.00	17.77	C
ATOM	2640	CG	ARG A 192	-53.822	30.683	6.881	1.00	18.47	C
ATOM	2643	CD	ARG A 192	-53.082	30.748	8.196	1.00	19.80	C
ATOM	2646	NE	ARG A 192	-52.174	31.900	8.224	1.00	20.33	N
ATOM	2648	CZ	ARG A 192	-51.176	32.048	9.077	1.00	20.10	C
ATOM	2649	NH1	ARG A 192	-50.928	31.113	9.972	1.00	21.72	N
ATOM	2652	NH2	ARG A 192	-50.407	33.119	9.019	1.00	20.79	N
ATOM	2655	C	ARG A 192	-56.477	28.311	5.590	1.00	18.00	C
ATOM	2656	O	ARG A 192	-56.067	27.170	5.871	1.00	18.21	O
ATOM	2658	N	ARG A 193	-57.758	28.588	5.335	1.00	17.59	N
ATOM	2659	CA	ARG A 193	-58.779	27.563	5.382	1.00	17.16	C
ATOM	2661	CB	ARG A 193	-60.141	28.164	5.741	1.00	17.45	C
ATOM	2664	CG	ARG A 193	-61.109	27.132	6.322	1.00	18.70	C
ATOM	2667	CD	ARG A 193	-62.479	27.682	6.686	1.00	20.20	C
ATOM	2670	NE	ARG A 193	-62.423	28.682	7.751	1.00	22.59	N
ATOM	2672	CZ	ARG A 193	-62.386	28.435	9.067	1.00	24.12	C
ATOM	2673	NH1	ARG A 193	-62.395	27.198	9.575	1.00	22.66	N
ATOM	2676	NH2	ARG A 193	-62.348	29.469	9.898	1.00	26.36	N
ATOM	2679	C	ARG A 193	-58.871	26.842	4.058	1.00	16.63	C
ATOM	2680	O	ARG A 193	-58.861	27.461	3.005	1.00	15.93	O
ATOM	2682	N	THR A 194	-58.955	25.518	4.149	1.00	16.72	N
ATOM	2683	CA	THR A 194	-59.332	24.638	3.051	1.00	16.88	C
ATOM	2685	CB	THR A 194	-59.520	23.250	3.592	1.00	16.43	C
ATOM	2687	OG1	THR A 194	-58.252	22.605	3.583	1.00	17.74	O
ATOM	2689	CG2	THR A 194	-60.454	22.441	2.754	1.00	18.43	C
ATOM	2693	C	THR A 194	-60.607	25.112	2.384	1.00	17.28	C
ATOM	2694	O	THR A 194	-61.441	25.734	3.022	1.00	17.44	O
ATOM	2696	N	GLN A 195	-60.765	24.855	1.091	1.00	18.05	N
ATOM	2697	CA	GLN A 195	-61.925	25.402	.397	1.00	18.76	C
ATOM	2699	CB	GLN A 195	-61.768	25.381	-1.114	1.00	18.87	C
ATOM	2702	CG	GLN A 195	-63.145	25.523	-1.761	1.00	20.75	C
ATOM	2705	CD	GLN A 195	-63.095	25.922	-3.179	1.00	23.41	C
ATOM	2706	OE1	GLN A 195	-62.031	26.088	-3.733	1.00	27.07	O

ATOM	2707	NE2	GLN A 195	-64.248	26.090	-3.791	1.00	25.02	N
ATOM	2710	C	GLN A 195	-63.243	24.704	.765	1.00	18.93	C
ATOM	2711	O	GLN A 195	-64.169	25.333	1.291	1.00	19.42	O
ATOM	2713	N	ARG A 196	-63.345	23.418	.443	1.00	18.90	N
ATOM	2714	CA	ARG A 196	-64.577	22.685	.663	1.00	18.68	C
ATOM	2716	CB	ARG A 196	-64.356	21.186	.454	1.00	18.53	C
ATOM	2719	CG	ARG A 196	-64.851	20.644	-.888	1.00	18.31	C
ATOM	2722	CD	ARG A 196	-64.705	21.622	-2.037	1.00	17.93	C
ATOM	2725	NE	ARG A 196	-65.930	21.720	-2.834	1.00	18.52	N
ATOM	2727	CZ	ARG A 196	-66.480	22.857	-3.256	1.00	18.63	C
ATOM	2728	NH1	ARG A 196	-65.938	24.014	-2.937	1.00	19.84	N
ATOM	2731	NH2	ARG A 196	-67.582	22.844	-3.999	1.00	18.53	N
ATOM	2734	C	ARG A 196	-65.070	22.982	2.061	1.00	19.02	C
ATOM	2735	O	ARG A 196	-66.276	23.085	2.284	1.00	19.30	O
ATOM	2737	N	LEU A 197	-64.127	23.152	2.990	1.00	19.09	N
ATOM	2738	CA	LEU A 197	-64.444	23.432	4.379	1.00	19.19	C
ATOM	2740	CB	LEU A 197	-63.203	23.346	5.262	1.00	19.14	C
ATOM	2743	CG	LEU A 197	-63.135	22.097	6.130	1.00	19.99	C
ATOM	2745	CD1	LEU A 197	-64.476	21.899	6.811	1.00	21.56	C
ATOM	2749	CD2	LEU A 197	-62.016	22.198	7.161	1.00	20.06	C
ATOM	2753	C	LEU A 197	-65.071	24.781	4.553	1.00	19.52	C
ATOM	2754	O	LEU A 197	-66.052	24.902	5.252	1.00	19.85	O
ATOM	2756	N	GLU A 198	-64.488	25.803	3.941	1.00	20.11	N
ATOM	2757	CA	GLU A 198	-65.051	27.147	3.986	1.00	20.47	C
ATOM	2759	CB	GLU A 198	-64.079	28.156	3.385	1.00	20.84	C
ATOM	2762	CG	GLU A 198	-64.659	29.539	3.069	1.00	22.30	C
ATOM	2765	CD	GLU A 198	-65.206	30.260	4.283	1.00	24.40	C
ATOM	2766	OE1	GLU A 198	-64.852	29.894	5.418	1.00	25.34	O
ATOM	2767	OE2	GLU A 198	-65.991	31.217	4.106	1.00	26.83	O
ATOM	2768	C	GLU A 198	-66.360	27.186	3.226	1.00	20.57	C
ATOM	2769	O	GLU A 198	-67.212	28.009	3.529	1.00	20.65	O
ATOM	2771	N	ALA A 199	-66.519	26.298	2.240	1.00	20.76	N
ATOM	2772	CA	ALA A 199	-67.751	26.229	1.452	1.00	20.52	C
ATOM	2774	CB	ALA A 199	-67.521	25.492	.166	1.00	20.44	C
ATOM	2778	C	ALA A 199	-68.907	25.611	2.232	1.00	20.55	C
ATOM	2779	O	ALA A 199	-69.988	26.158	2.202	1.00	20.17	O
ATOM	2781	N	VAL A 200	-68.701	24.493	2.933	1.00	21.10	N
ATOM	2782	CA	VAL A 200	-69.791	23.935	3.774	1.00	21.51	C
ATOM	2784	CB	VAL A 200	-69.509	22.576	4.485	1.00	21.18	C
ATOM	2786	CG1	VAL A 200	-69.618	21.460	3.525	1.00	21.27	C
ATOM	2790	CG2	VAL A 200	-68.161	22.569	5.200	1.00	21.23	C
ATOM	2794	C	VAL A 200	-70.170	24.905	4.860	1.00	21.99	C
ATOM	2795	O	VAL A 200	-71.338	25.022	5.186	1.00	22.99	O
ATOM	2797	N	TRP A 201	-69.204	25.590	5.440	1.00	22.10	N
ATOM	2798	CA	TRP A 201	-69.551	26.565	6.433	1.00	22.58	C
ATOM	2800	CB	TRP A 201	-68.315	27.144	7.121	1.00	22.92	C
ATOM	2803	CG	TRP A 201	-68.667	27.957	8.310	1.00	23.65	C
ATOM	2804	CD1	TRP A 201	-68.758	27.520	9.588	1.00	25.08	C
ATOM	2806	NE1	TRP A 201	-69.127	28.549	10.417	1.00	26.35	N
ATOM	2808	CE2	TRP A 201	-69.290	29.683	9.669	1.00	25.48	C
ATOM	2809	CD2	TRP A 201	-69.010	29.344	8.332	1.00	24.77	C
ATOM	2810	CE3	TRP A 201	-69.095	30.336	7.355	1.00	24.92	C
ATOM	2812	CZ3	TRP A 201	-69.450	31.616	7.739	1.00	24.96	C
ATOM	2814	CH2	TRP A 201	-69.720	31.918	9.078	1.00	25.08	C
ATOM	2816	CZ2	TRP A 201	-69.648	30.968	10.056	1.00	24.68	C
ATOM	2818	C	TRP A 201	-70.351	27.662	5.758	1.00	22.71	C
ATOM	2819	O	TRP A 201	-71.482	27.921	6.140	1.00	22.76	O
ATOM	2821	N	SER A 202	-69.776	28.272	4.727	1.00	23.04	N
ATOM	2822	CA	SER A 202	-70.335	29.500	4.165	1.00	23.29	C
ATOM	2824	CB	SER A 202	-69.304	30.205	3.290	1.00	23.20	C

ATOM	2827	OG SER A 202	-68.306	30.807	4.111	1.00	23.07	O
ATOM	2829	C SER A 202	-71.666	29.341	3.432	1.00	23.66	C
ATOM	2830	O SER A 202	-72.392	30.310	3.265	1.00	23.38	O
ATOM	2832	N ILE A 203	-71.996	28.124	3.016	1.00	24.49	N
ATOM	2833	CA ILE A 203	-73.321	27.859	2.460	1.00	25.07	C
ATOM	2835	CB ILE A 203	-73.391	26.530	1.653	1.00	24.74	C
ATOM	2837	CG1 ILE A 203	-72.494	26.587	.423	1.00	24.36	C
ATOM	2840	CD1 ILE A 203	-72.344	25.255	-.281	1.00	23.81	C
ATOM	2844	CG2 ILE A 203	-74.795	26.292	1.156	1.00	24.45	C
ATOM	2848	C ILE A 203	-74.338	27.858	3.607	1.00	25.91	C
ATOM	2849	O ILE A 203	-75.350	28.554	3.546	1.00	26.26	O
ATOM	2851	N GLU A 204	-74.051	27.104	4.662	1.00	26.69	N
ATOM	2852	CA GLU A 204	-74.929	27.053	5.824	1.00	27.40	C
ATOM	2854	CB GLU A 204	-74.286	26.217	6.931	1.00	27.65	C
ATOM	2857	CG GLU A 204	-75.132	26.032	8.183	1.00	29.02	C
ATOM	2860	CD GLU A 204	-76.530	25.491	7.903	1.00	30.44	C
ATOM	2861	OE1 GLU A 204	-76.746	24.858	6.844	1.00	31.22	O
ATOM	2862	OE2 GLU A 204	-77.417	25.703	8.753	1.00	31.96	O
ATOM	2863	C GLU A 204	-75.258	28.451	6.342	1.00	27.61	C
ATOM	2864	O GLU A 204	-76.409	28.778	6.526	1.00	27.78	O
ATOM	2866	N ALA A 205	-74.240	29.274	6.547	1.00	28.15	N
ATOM	2867	CA ALA A 205	-74.416	30.637	7.050	1.00	28.41	C
ATOM	2869	CB ALA A 205	-73.062	31.285	7.274	1.00	28.41	C
ATOM	2873	C ALA A 205	-75.234	31.520	6.136	1.00	28.84	C
ATOM	2874	O ALA A 205	-76.050	32.302	6.610	1.00	28.54	O
ATOM	2876	N TYR A 206	-74.978	31.408	4.832	1.00	29.84	N
ATOM	2877	CA TYR A 206	-75.581	32.281	3.805	1.00	30.41	C
ATOM	2879	CB TYR A 206	-74.801	32.170	2.493	1.00	30.20	C
ATOM	2882	CG TYR A 206	-75.180	33.173	1.425	1.00	29.62	C
ATOM	2883	CD1 TYR A 206	-74.721	34.482	1.484	1.00	30.32	C
ATOM	2885	CE1 TYR A 206	-75.047	35.412	.489	1.00	29.70	C
ATOM	2887	CZ TYR A 206	-75.830	35.028	-.578	1.00	28.41	C
ATOM	2888	OH TYR A 206	-76.155	35.956	-1.548	1.00	26.74	O
ATOM	2890	CE2 TYR A 206	-76.291	33.724	-.659	1.00	28.16	C
ATOM	2892	CD2 TYR A 206	-75.958	32.806	.335	1.00	28.39	C
ATOM	2894	C TYR A 206	-77.038	31.918	3.561	1.00	31.15	C
ATOM	2895	O TYR A 206	-77.858	32.781	3.268	1.00	31.30	O
ATOM	2897	N ARG A 207	-77.341	30.631	3.676	1.00	32.05	N
ATOM	2898	CA ARG A 207	-78.691	30.128	3.555	1.00	32.81	C
ATOM	2900	CB ARG A 207	-78.684	28.623	3.813	1.00	32.86	C
ATOM	2903	CG ARG A 207	-80.028	27.941	3.673	1.00	34.01	C
ATOM	2906	CD ARG A 207	-80.066	26.632	4.450	1.00	35.11	C
ATOM	2909	NE ARG A 207	-79.096	25.669	3.939	1.00	35.35	N
ATOM	2911	CZ ARG A 207	-79.353	24.737	3.022	1.00	36.00	C
ATOM	2912	NH1 ARG A 207	-80.562	24.611	2.478	1.00	35.55	N
ATOM	2915	NH2 ARG A 207	-78.381	23.918	2.640	1.00	37.05	N
ATOM	2918	C ARG A 207	-79.609	30.826	4.550	1.00	33.57	C
ATOM	2919	O ARG A 207	-80.783	31.003	4.263	1.00	33.82	O
ATOM	2921	N LYS A 208	-79.064	31.223	5.707	1.00	34.58	N
ATOM	2922	CA LYS A 208	-79.835	31.830	6.812	1.00	35.21	C
ATOM	2924	CB LYS A 208	-79.137	31.599	8.159	1.00	35.05	C
ATOM	2927	CG LYS A 208	-78.847	30.147	8.496	1.00	34.78	C
ATOM	2930	CD LYS A 208	-78.583	29.974	9.980	1.00	34.44	C
ATOM	2933	CE LYS A 208	-77.992	28.619	10.309	1.00	34.16	C
ATOM	2936	NZ LYS A 208	-76.540	28.693	10.572	1.00	33.81	N
ATOM	2940	C LYS A 208	-80.075	33.333	6.663	1.00	36.06	C
ATOM	2941	O LYS A 208	-81.032	33.857	7.226	1.00	36.19	O
ATOM	2943	N LYS A 209	-79.197	34.032	5.946	1.00	36.93	N
ATOM	2944	CA LYS A 209	-79.412	35.445	5.652	1.00	37.77	C
ATOM	2946	CB LYS A 209	-78.284	35.998	4.779	1.00	38.18	C

ATOM	2949	CG	LYS	A 209	-76.928	36.134	5.458	1.00	39.90	C
ATOM	2952	CD	LYS	A 209	-75.982	37.023	4.621	1.00	42.39	C
ATOM	2955	CE	LYS	A 209	-74.531	37.016	5.169	1.00	44.21	C
ATOM	2958	NZ	LYS	A 209	-74.350	37.688	6.509	1.00	44.75	N
ATOM	2962	C	LYS	A 209	-80.738	35.600	4.911	1.00	37.89	C
ATOM	2963	O	LYS	A 209	-81.049	34.785	4.036	1.00	38.20	O
ATOM	2965	N	GLU	A 210	-81.518	36.632	5.243	1.00	37.87	N
ATOM	2966	CA	GLU	A 210	-82.847	36.797	4.630	1.00	37.72	C
ATOM	2968	CB	GLU	A 210	-83.821	37.578	5.549	1.00	38.19	C
ATOM	2971	CG	GLU	A 210	-83.850	39.130	5.404	1.00	39.68	C
ATOM	2974	CD	GLU	A 210	-85.270	39.736	5.534	1.00	41.27	C
ATOM	2975	OE1	GLU	A 210	-86.172	39.100	6.136	1.00	41.67	O
ATOM	2976	OE2	GLU	A 210	-85.482	40.859	5.019	1.00	42.17	O
ATOM	2977	C	GLU	A 210	-82.760	37.394	3.219	1.00	36.68	C
ATOM	2978	O	GLU	A 210	-83.605	37.103	2.371	1.00	36.39	O
ATOM	2980	N	ASP	A 211	-81.725	38.198	2.969	1.00	35.65	N
ATOM	2981	CA	ASP	A 211	-81.481	38.777	1.632	1.00	34.94	C
ATOM	2983	CB	ASP	A 211	-81.064	40.257	1.740	1.00	35.22	C
ATOM	2986	CG	ASP	A 211	-79.811	40.477	2.599	1.00	35.84	C
ATOM	2987	OD1	ASP	A 211	-79.379	39.557	3.341	1.00	35.60	O
ATOM	2988	OD2	ASP	A 211	-79.276	41.606	2.536	1.00	36.73	O
ATOM	2989	C	ASP	A 211	-80.451	37.989	.811	1.00	33.73	C
ATOM	2990	O	ASP	A 211	-79.753	38.556	-.021	1.00	33.40	O
ATOM	2992	N	ALA	A 212	-80.366	36.682	1.057	1.00	32.53	N
ATOM	2993	CA	ALA	A 212	-79.506	35.783	.293	1.00	31.40	C
ATOM	2995	CB	ALA	A 212	-79.569	34.382	.889	1.00	30.79	C
ATOM	2999	C	ALA	A 212	-79.940	35.762	-1.176	1.00	30.60	C
ATOM	3000	O	ALA	A 212	-81.124	35.594	-1.467	1.00	30.37	O
ATOM	3002	N	ASN	A 213	-79.000	35.967	-2.097	1.00	29.71	N
ATOM	3003	CA	ASN	A 213	-79.304	35.809	-3.520	1.00	29.27	C
ATOM	3005	CB	ASN	A 213	-78.144	36.315	-4.402	1.00	29.30	C
ATOM	3008	CG	ASN	A 213	-78.442	36.217	-5.911	1.00	29.07	C
ATOM	3009	OD1	ASN	A 213	-78.437	37.215	-6.624	1.00	29.38	O
ATOM	3010	ND2	ASN	A 213	-78.691	35.015	-6.389	1.00	28.37	N
ATOM	3013	C	ASN	A 213	-79.634	34.328	-3.785	1.00	28.81	C
ATOM	3014	O	ASN	A 213	-78.784	33.455	-3.658	1.00	29.04	O
ATOM	3016	N	GLN	A 214	-80.880	34.053	-4.145	1.00	28.04	N
ATOM	3017	CA	GLN	A 214	-81.376	32.692	-4.191	1.00	27.32	C
ATOM	3019	CB	GLN	A 214	-82.906	32.700	-4.178	1.00	27.47	C
ATOM	3022	CG	GLN	A 214	-83.505	33.269	-2.889	1.00	27.68	C
ATOM	3025	CD	GLN	A 214	-83.294	32.359	-1.688	1.00	27.54	C
ATOM	3026	OE1	GLN	A 214	-83.940	31.320	-1.571	1.00	28.19	O
ATOM	3027	NE2	GLN	A 214	-82.392	32.748	-.790	1.00	26.90	N
ATOM	3030	C	GLN	A 214	-80.858	31.902	-5.385	1.00	26.71	C
ATOM	3031	O	GLN	A 214	-80.863	30.668	-5.355	1.00	26.91	O
ATOM	3033	N	VAL	A 215	-80.439	32.600	-6.442	1.00	25.81	N
ATOM	3034	CA	VAL	A 215	-79.789	31.947	-7.588	1.00	24.88	C
ATOM	3036	CB	VAL	A 215	-79.651	32.894	-8.808	1.00	24.81	C
ATOM	3038	CG1	VAL	A 215	-78.459	32.483	-9.680	1.00	24.94	C
ATOM	3042	CG2	VAL	A 215	-80.935	32.907	-9.620	1.00	23.74	C
ATOM	3046	C	VAL	A 215	-78.415	31.427	-7.176	1.00	24.18	C
ATOM	3047	O	VAL	A 215	-78.073	30.282	-7.430	1.00	23.74	O
ATOM	3049	N	LEU	A 216	-77.650	32.283	-6.510	1.00	23.78	N
ATOM	3050	CA	LEU	A 216	-76.299	31.961	-6.036	1.00	23.33	C
ATOM	3052	CB	LEU	A 216	-75.680	33.195	-5.373	1.00	23.10	C
ATOM	3055	CG	LEU	A 216	-74.194	33.192	-5.059	1.00	22.91	C
ATOM	3057	CD1	LEU	A 216	-73.383	32.962	-6.312	1.00	22.88	C
ATOM	3061	CD2	LEU	A 216	-73.826	34.518	-4.424	1.00	22.03	C
ATOM	3065	C	LEU	A 216	-76.333	30.797	-5.050	1.00	22.92	C
ATOM	3066	O	LEU	A 216	-75.598	29.826	-5.201	1.00	23.38	O

ATOM	3068	N	LEU A 217	-77.214	30.890	-4.060	1.00	22.09	N
ATOM	3069	CA	LEU A 217	-77.367	29.851	-3.054	1.00	21.14	C
ATOM	3071	CB	LEU A 217	-78.488	30.225	-2.090	1.00	21.11	C
ATOM	3074	CG	LEU A 217	-78.782	29.236	-.961	1.00	21.08	C
ATOM	3076	CD1	LEU A 217	-77.703	29.267	.101	1.00	19.87	C
ATOM	3080	CD2	LEU A 217	-80.146	29.535	-.360	1.00	21.46	C
ATOM	3084	C	LEU A 217	-77.682	28.517	-3.685	1.00	20.38	C
ATOM	3085	O	LEU A 217	-77.133	27.508	-3.286	1.00	20.04	O
ATOM	3087	N	GLU A 218	-78.579	28.512	-4.661	1.00	19.97	N
ATOM	3088	CA	GLU A 218	-79.027	27.260	-5.273	1.00	19.72	C
ATOM	3090	CB	GLU A 218	-80.220	27.499	-6.212	1.00	19.87	C
ATOM	3093	CG	GLU A 218	-81.005	26.228	-6.627	1.00	21.03	C
ATOM	3096	CD	GLU A 218	-82.303	26.535	-7.424	1.00	22.59	C
ATOM	3097	OE1	GLU A 218	-82.643	27.730	-7.619	1.00	23.61	O
ATOM	3098	OE2	GLU A 218	-82.983	25.578	-7.863	1.00	22.31	O
ATOM	3099	C	GLU A 218	-77.854	26.613	-6.009	1.00	19.03	C
ATOM	3100	O	GLU A 218	-77.522	25.452	-5.764	1.00	18.97	O
ATOM	3102	N	LEU A 219	-77.203	27.381	-6.878	1.00	18.06	N
ATOM	3103	CA	LEU A 219	-76.012	26.907	-7.573	1.00	17.15	C
ATOM	3105	CB	LEU A 219	-75.439	28.029	-8.418	1.00	16.75	C
ATOM	3108	CG	LEU A 219	-74.196	27.717	-9.231	1.00	16.01	C
ATOM	3110	CD1	LEU A 219	-74.511	26.832	-10.404	1.00	12.99	C
ATOM	3114	CD2	LEU A 219	-73.577	29.028	-9.675	1.00	16.20	C
ATOM	3118	C	LEU A 219	-74.973	26.409	-6.570	1.00	16.75	C
ATOM	3119	O	LEU A 219	-74.377	25.355	-6.757	1.00	16.47	O
ATOM	3121	N	ALA A 220	-74.786	27.169	-5.496	1.00	16.40	N
ATOM	3122	CA	ALA A 220	-73.858	26.808	-4.422	1.00	16.15	C
ATOM	3124	CB	ALA A 220	-73.906	27.830	-3.302	1.00	16.13	C
ATOM	3128	C	ALA A 220	-74.137	25.436	-3.862	1.00	16.03	C
ATOM	3129	O	ALA A 220	-73.253	24.603	-3.799	1.00	16.30	O
ATOM	3131	N	ILE A 221	-75.371	25.201	-3.447	1.00	16.24	N
ATOM	3132	CA	ILE A 221	-75.731	23.916	-2.857	1.00	16.33	C
ATOM	3134	CB	ILE A 221	-77.203	23.895	-2.352	1.00	15.94	C
ATOM	3136	CG1	ILE A 221	-77.410	24.879	-1.203	1.00	15.36	C
ATOM	3139	CD1	ILE A 221	-78.842	25.234	-.962	1.00	14.49	C
ATOM	3143	CG2	ILE A 221	-77.571	22.537	-1.836	1.00	15.18	C
ATOM	3147	C	ILE A 221	-75.509	22.841	-3.911	1.00	17.06	C
ATOM	3148	O	ILE A 221	-74.899	21.819	-3.653	1.00	16.51	O
ATOM	3150	N	LEU A 222	-75.981	23.138	-5.117	1.00	18.53	N
ATOM	3151	CA	LEU A 222	-76.006	22.207	-6.243	1.00	19.25	C
ATOM	3153	CB	LEU A 222	-76.657	22.891	-7.445	1.00	19.18	C
ATOM	3156	CG	LEU A 222	-77.037	21.972	-8.597	1.00	19.65	C
ATOM	3158	CD1	LEU A 222	-78.318	22.489	-9.268	1.00	19.24	C
ATOM	3162	CD2	LEU A 222	-75.869	21.788	-9.597	1.00	18.75	C
ATOM	3166	C	LEU A 222	-74.618	21.751	-6.630	1.00	20.07	C
ATOM	3167	O	LEU A 222	-74.381	20.558	-6.778	1.00	20.55	O
ATOM	3169	N	ASP A 223	-73.710	22.711	-6.797	1.00	20.92	N
ATOM	3170	CA	ASP A 223	-72.330	22.427	-7.189	1.00	21.43	C
ATOM	3172	CB	ASP A 223	-71.621	23.723	-7.585	1.00	21.63	C
ATOM	3175	CG	ASP A 223	-70.202	23.495	-8.070	1.00	22.51	C
ATOM	3176	OD1	ASP A 223	-69.296	23.282	-7.225	1.00	24.67	O
ATOM	3177	OD2	ASP A 223	-69.990	23.544	-9.296	1.00	22.73	O
ATOM	3178	C	ASP A 223	-71.573	21.722	-6.060	1.00	21.79	C
ATOM	3179	O	ASP A 223	-70.843	20.764	-6.313	1.00	21.73	O
ATOM	3181	N	TYR A 224	-71.755	22.182	-4.820	1.00	22.14	N
ATOM	3182	CA	TYR A 224	-71.120	21.529	-3.686	1.00	22.22	C
ATOM	3184	CB	TYR A 224	-71.462	22.190	-2.355	1.00	22.38	C
ATOM	3187	CG	TYR A 224	-70.673	21.542	-1.249	1.00	22.65	C
ATOM	3188	CD1	TYR A 224	-69.356	21.864	-1.053	1.00	22.59	C
ATOM	3190	CE1	TYR A 224	-68.622	21.258	-.098	1.00	23.38	C

ATOM	3192	CZ	TYR A 224	-69.177	20.288	.670	1.00	23.22	C
ATOM	3193	OH	TYR A 224	-68.402	19.684	1.620	1.00	23.93	O
ATOM	3195	CE2	TYR A 224	-70.483	19.931	.500	1.00	23.36	C
ATOM	3197	CD2	TYR A 224	-71.221	20.547	-.464	1.00	23.41	C
ATOM	3199	C	TYR A 224	-71.485	20.065	-3.605	1.00	22.52	C
ATOM	3200	O	TYR A 224	-70.641	19.233	-3.311	1.00	22.93	O
ATOM	3202	N	ASN A 225	-72.742	19.749	-3.856	1.00	22.94	N
ATOM	3203	CA	ASN A 225	-73.187	18.374	-3.801	1.00	23.44	C
ATOM	3205	CB	ASN A 225	-74.706	18.310	-3.729	1.00	23.34	C
ATOM	3208	CG	ASN A 225	-75.210	18.579	-2.352	1.00	22.62	C
ATOM	3209	OD1	ASN A 225	-74.749	17.977	-1.407	1.00	23.33	O
ATOM	3210	ND2	ASN A 225	-76.146	19.494	-2.223	1.00	22.52	N
ATOM	3213	C	ASN A 225	-72.694	17.523	-4.959	1.00	24.26	C
ATOM	3214	O	ASN A 225	-72.408	16.328	-4.767	1.00	24.88	O
ATOM	3216	N	MET A 226	-72.607	18.119	-6.150	1.00	24.67	N
ATOM	3217	CA	MET A 226	-72.167	17.395	-7.349	1.00	24.95	C
ATOM	3219	CB	MET A 226	-72.421	18.237	-8.594	1.00	25.37	C
ATOM	3222	CG	MET A 226	-71.785	17.713	-9.873	1.00	27.23	C
ATOM	3225	SD	MET A 226	-70.837	19.010	-10.708	1.00	31.73	S
ATOM	3226	CE	MET A 226	-72.185	20.000	-11.367	1.00	31.30	C
ATOM	3230	C	MET A 226	-70.688	17.045	-7.247	1.00	24.70	C
ATOM	3231	O	MET A 226	-70.291	15.929	-7.568	1.00	24.91	O
ATOM	3233	N	ILE A 227	-69.873	17.995	-6.799	1.00	24.43	N
ATOM	3234	CA	ILE A 227	-68.456	17.726	-6.598	1.00	24.18	C
ATOM	3236	CB	ILE A 227	-67.656	18.971	-6.154	1.00	24.15	C
ATOM	3238	CG1	ILE A 227	-67.628	20.021	-7.253	1.00	23.64	C
ATOM	3241	CD1	ILE A 227	-66.853	21.246	-6.880	1.00	23.11	C
ATOM	3245	CG2	ILE A 227	-66.229	18.603	-5.827	1.00	24.36	C
ATOM	3249	C	ILE A 227	-68.329	16.645	-5.546	1.00	23.98	C
ATOM	3250	O	ILE A 227	-67.609	15.689	-5.747	1.00	24.35	O
ATOM	3252	N	GLN A 228	-69.048	16.779	-4.439	1.00	23.79	N
ATOM	3253	CA	GLN A 228	-69.020	15.754	-3.382	1.00	23.66	C
ATOM	3255	CB	GLN A 228	-70.077	16.010	-2.305	1.00	23.64	C
ATOM	3258	CG	GLN A 228	-69.972	15.049	-1.145	1.00	22.91	C
ATOM	3261	CD	GLN A 228	-70.891	15.417	-.033	1.00	23.26	C
ATOM	3262	OE1	GLN A 228	-72.082	15.098	-.069	1.00	25.64	O
ATOM	3263	NE2	GLN A 228	-70.357	16.083	.978	1.00	21.77	N
ATOM	3266	C	GLN A 228	-69.243	14.349	-3.893	1.00	23.54	C
ATOM	3267	O	GLN A 228	-68.664	13.412	-3.357	1.00	23.59	O
ATOM	3269	N	SER A 229	-70.105	14.194	-4.893	1.00	23.27	N
ATOM	3270	CA	SER A 229	-70.394	12.868	-5.406	1.00	23.35	C
ATOM	3272	CB	SER A 229	-71.753	12.812	-6.095	1.00	23.31	C
ATOM	3275	OG	SER A 229	-71.836	13.823	-7.060	1.00	24.31	O
ATOM	3277	C	SER A 229	-69.287	12.394	-6.331	1.00	23.19	C
ATOM	3278	O	SER A 229	-69.130	11.194	-6.512	1.00	23.44	O
ATOM	3280	N	VAL A 230	-68.512	13.306	-6.914	1.00	23.14	N
ATOM	3281	CA	VAL A 230	-67.269	12.874	-7.552	1.00	23.00	C
ATOM	3283	CB	VAL A 230	-66.469	13.998	-8.253	1.00	22.75	C
ATOM	3285	CG1	VAL A 230	-65.091	13.470	-8.667	1.00	21.87	C
ATOM	3289	CG2	VAL A 230	-67.222	14.539	-9.459	1.00	21.74	C
ATOM	3293	C	VAL A 230	-66.417	12.268	-6.452	1.00	23.35	C
ATOM	3294	O	VAL A 230	-65.917	11.173	-6.589	1.00	23.68	O
ATOM	3296	N	TYR A 231	-66.284	12.968	-5.340	1.00	23.96	N
ATOM	3297	CA	TYR A 231	-65.414	12.509	-4.274	1.00	24.60	C
ATOM	3299	CB	TYR A 231	-65.486	13.426	-3.051	1.00	24.69	C
ATOM	3302	CG	TYR A 231	-64.963	14.837	-3.211	1.00	24.11	C
ATOM	3303	CD1	TYR A 231	-64.071	15.190	-4.214	1.00	24.31	C
ATOM	3305	CE1	TYR A 231	-63.592	16.491	-4.322	1.00	24.63	C
ATOM	3307	CZ	TYR A 231	-63.994	17.441	-3.406	1.00	25.13	C
ATOM	3308	OH	TYR A 231	-63.554	18.754	-3.456	1.00	24.64	O

ATOM	3310	CE2 TYR A 231	-64.863	17.086	-2.402	1.00	25.76	C
ATOM	3312	CD2 TYR A 231	-65.329	15.800	-2.308	1.00	24.28	C
ATOM	3314	C TYR A 231	-65.785	11.111	-3.833	1.00	25.26	C
ATOM	3315	O TYR A 231	-64.917	10.315	-3.488	1.00	25.44	O
ATOM	3317	N GLN A 232	-67.079	10.819	-3.833	1.00	26.04	N
ATOM	3318	CA GLN A 232	-67.566	9.532	-3.362	1.00	26.57	C
ATOM	3320	CB GLN A 232	-69.060	9.621	-3.057	1.00	26.44	C
ATOM	3323	CG GLN A 232	-69.339	10.397	-1.778	1.00	26.25	C
ATOM	3326	CD GLN A 232	-70.786	10.819	-1.630	1.00	26.53	C
ATOM	3327	OE1 GLN A 232	-71.667	10.305	-2.318	1.00	27.77	O
ATOM	3328	NE2 GLN A 232	-71.040	11.759	-.723	1.00	25.51	N
ATOM	3331	C GLN A 232	-67.238	8.414	-4.351	1.00	27.44	C
ATOM	3332	O GLN A 232	-66.852	7.319	-3.935	1.00	27.33	O
ATOM	3334	N ARG A 233	-67.373	8.691	-5.649	1.00	28.66	N
ATOM	3335	CA ARG A 233	-66.905	7.768	-6.685	1.00	29.86	C
ATOM	3337	CB ARG A 233	-67.212	8.292	-8.090	1.00	30.08	C
ATOM	3340	CG ARG A 233	-66.378	7.634	-9.179	1.00	32.80	C
ATOM	3343	CD ARG A 233	-66.913	7.923	-10.572	1.00	36.74	C
ATOM	3346	NE ARG A 233	-66.962	9.363	-10.860	1.00	40.33	N
ATOM	3348	CZ ARG A 233	-68.069	10.119	-10.891	1.00	43.12	C
ATOM	3349	NH1 ARG A 233	-69.276	9.600	-10.655	1.00	44.55	N
ATOM	3352	NH2 ARG A 233	-67.971	11.418	-11.169	1.00	43.38	N
ATOM	3355	C ARG A 233	-65.401	7.518	-6.514	1.00	30.34	C
ATOM	3356	O ARG A 233	-64.962	6.374	-6.412	1.00	30.55	O
ATOM	3358	N ASP A 234	-64.621	8.591	-6.461	1.00	31.00	N
ATOM	3359	CA ASP A 234	-63.187	8.487	-6.205	1.00	31.48	C
ATOM	3361	CB ASP A 234	-62.594	9.870	-5.942	1.00	31.56	C
ATOM	3364	CG ASP A 234	-62.573	10.741	-7.167	1.00	32.78	C
ATOM	3365	OD1 ASP A 234	-62.911	10.249	-8.276	1.00	34.25	O
ATOM	3366	OD2 ASP A 234	-62.215	11.930	-7.011	1.00	34.81	O
ATOM	3367	C ASP A 234	-62.885	7.596	-5.010	1.00	31.70	C
ATOM	3368	O ASP A 234	-61.985	6.780	-5.053	1.00	31.56	O
ATOM	3370	N LEU A 235	-63.647	7.766	-3.943	1.00	32.43	N
ATOM	3371	CA LEU A 235	-63.366	7.102	-2.685	1.00	33.04	C
ATOM	3373	CB LEU A 235	-64.072	7.833	-1.553	1.00	32.57	C
ATOM	3376	CG LEU A 235	-63.884	7.252	-.167	1.00	30.85	C
ATOM	3378	CD1 LEU A 235	-62.428	7.261	.174	1.00	28.78	C
ATOM	3382	CD2 LEU A 235	-64.703	8.050	.823	1.00	29.93	C
ATOM	3386	C LEU A 235	-63.789	5.637	-2.686	1.00	34.69	C
ATOM	3387	O LEU A 235	-63.123	4.812	-2.065	1.00	35.29	O
ATOM	3389	N ARG A 236	-64.893	5.299	-3.354	1.00	36.09	N
ATOM	3390	CA ARG A 236	-65.290	3.893	-3.458	1.00	37.26	C
ATOM	3392	CB ARG A 236	-66.672	3.736	-4.103	1.00	37.64	C
ATOM	3395	CG ARG A 236	-67.839	4.000	-3.157	1.00	39.15	C
ATOM	3398	CD ARG A 236	-69.162	3.528	-3.741	1.00	40.25	C
ATOM	3401	NE ARG A 236	-69.385	4.046	-5.095	1.00	41.64	N
ATOM	3403	CZ ARG A 236	-69.837	5.269	-5.392	1.00	42.54	C
ATOM	3404	NH1 ARG A 236	-70.115	6.147	-4.430	1.00	42.18	N
ATOM	3407	NH2 ARG A 236	-70.012	5.621	-6.668	1.00	43.04	N
ATOM	3410	C ARG A 236	-64.266	3.108	-4.262	1.00	37.81	C
ATOM	3411	O ARG A 236	-64.104	1.912	-4.060	1.00	38.07	O
ATOM	3413	N GLU A 237	-63.584	3.801	-5.168	1.00	38.69	N
ATOM	3414	CA GLU A 237	-62.639	3.189	-6.102	1.00	39.45	C
ATOM	3416	CB GLU A 237	-62.452	4.108	-7.338	1.00	40.25	C
ATOM	3419	CG GLU A 237	-62.268	3.387	-8.710	1.00	43.02	C
ATOM	3422	CD GLU A 237	-60.835	3.458	-9.287	1.00	46.64	C
ATOM	3423	OE1 GLU A 237	-60.093	4.440	-8.997	1.00	48.03	O
ATOM	3424	OE2 GLU A 237	-60.466	2.519	-10.047	1.00	48.35	O
ATOM	3425	C GLU A 237	-61.310	2.905	-5.408	1.00	38.70	C
ATOM	3426	O GLU A 237	-60.832	1.782	-5.440	1.00	38.33	O

ATOM	3428	N	THR A 238	-60.726	3.908	-4.764	1.00	38.40	N
ATOM	3429	CA	THR A 238	-59.504	3.669	-4.009	1.00	38.74	C
ATOM	3431	CB	THR A 238	-58.710	4.965	-3.593	1.00	38.78	C
ATOM	3433	OG1	THR A 238	-59.294	5.583	-2.444	1.00	38.38	O
ATOM	3435	CG2	THR A 238	-58.612	5.963	-4.743	1.00	38.94	C
ATOM	3439	C	THR A 238	-59.802	2.834	-2.773	1.00	39.05	C
ATOM	3440	O	THR A 238	-58.886	2.294	-2.153	1.00	39.17	O
ATOM	3442	N	SER A 239	-61.076	2.717	-2.408	1.00	39.34	N
ATOM	3443	CA	SER A 239	-61.443	1.841	-1.304	1.00	39.47	C
ATOM	3445	CB	SER A 239	-62.844	2.149	-.792	1.00	39.34	C
ATOM	3448	OG	SER A 239	-63.087	1.428	.396	1.00	39.94	O
ATOM	3450	C	SER A 239	-61.313	.367	-1.699	1.00	39.53	C
ATOM	3451	O	SER A 239	-60.834	-.439	-.900	1.00	39.53	O
ATOM	3453	N	ARG A 240	-61.728	.021	-2.921	1.00	39.69	N
ATOM	3454	CA	ARG A 240	-61.559	-1.341	-3.439	1.00	39.79	C
ATOM	3456	CB	ARG A 240	-62.105	-1.492	-4.867	1.00	40.34	C
ATOM	3459	CG	ARG A 240	-63.624	-1.715	-4.947	1.00	43.03	C
ATOM	3462	CD	ARG A 240	-64.104	-2.125	-6.364	1.00	46.66	C
ATOM	3465	NE	ARG A 240	-63.780	-1.136	-7.416	1.00	50.46	N
ATOM	3467	CZ	ARG A 240	-64.535	-.082	-7.775	1.00	53.13	C
ATOM	3468	NH1	ARG A 240	-64.112	.732	-8.749	1.00	53.13	N
ATOM	3471	NH2	ARG A 240	-65.707	.178	-7.178	1.00	54.02	N
ATOM	3474	C	ARG A 240	-60.089	-1.704	-3.414	1.00	38.94	C
ATOM	3475	O	ARG A 240	-59.732	-2.795	-2.978	1.00	39.09	O
ATOM	3477	N	TRP A 241	-59.246	-.775	-3.862	1.00	37.92	N
ATOM	3478	CA	TRP A 241	-57.787	-.947	-3.842	1.00	36.93	C
ATOM	3480	CB	TRP A 241	-57.099	.318	-4.388	1.00	36.57	C
ATOM	3483	CG	TRP A 241	-55.624	.333	-4.209	1.00	34.71	C
ATOM	3484	CD1	TRP A 241	-54.702	-.312	-4.969	1.00	33.50	C
ATOM	3486	NE1	TRP A 241	-53.441	-.068	-4.482	1.00	32.62	N
ATOM	3488	CE2	TRP A 241	-53.538	.747	-3.392	1.00	31.16	C
ATOM	3489	CD2	TRP A 241	-54.896	1.022	-3.191	1.00	32.61	C
ATOM	3490	CE3	TRP A 241	-55.268	1.837	-2.117	1.00	32.63	C
ATOM	3492	CZ3	TRP A 241	-54.293	2.345	-1.307	1.00	31.31	C
ATOM	3494	CH2	TRP A 241	-52.956	2.053	-1.537	1.00	31.81	C
ATOM	3496	CZ2	TRP A 241	-52.560	1.252	-2.577	1.00	31.18	C
ATOM	3498	C	TRP A 241	-57.263	-1.274	-2.444	1.00	36.67	C
ATOM	3499	O	TRP A 241	-56.540	-2.247	-2.252	1.00	36.32	O
ATOM	3501	N	TRP A 242	-57.647	-.460	-1.471	1.00	36.59	N
ATOM	3502	CA	TRP A 242	-57.150	-.593	-.100	1.00	36.80	C
ATOM	3504	CB	TRP A 242	-57.694	.582	.734	1.00	36.47	C
ATOM	3507	CG	TRP A 242	-57.113	.754	2.107	1.00	35.59	C
ATOM	3508	CD1	TRP A 242	-57.806	.831	3.270	1.00	35.14	C
ATOM	3510	NE1	TRP A 242	-56.950	.999	4.328	1.00	34.47	N
ATOM	3512	CE2	TRP A 242	-55.667	1.032	3.859	1.00	34.64	C
ATOM	3513	CD2	TRP A 242	-55.728	.884	2.461	1.00	35.26	C
ATOM	3514	CE3	TRP A 242	-54.531	.882	1.732	1.00	35.22	C
ATOM	3516	CZ3	TRP A 242	-53.342	1.024	2.409	1.00	34.30	C
ATOM	3518	CH2	TRP A 242	-53.318	1.169	3.799	1.00	34.70	C
ATOM	3520	CZ2	TRP A 242	-54.467	1.173	4.542	1.00	34.63	C
ATOM	3522	C	TRP A 242	-57.482	-1.975	.524	1.00	37.34	C
ATOM	3523	O	TRP A 242	-56.628	-2.623	1.126	1.00	36.53	O
ATOM	3525	N	ARG A 243	-58.720	-2.421	.348	1.00	38.54	N
ATOM	3526	CA	ARG A 243	-59.149	-3.734	.822	1.00	39.71	C
ATOM	3528	CB	ARG A 243	-60.669	-3.896	.687	1.00	40.10	C
ATOM	3531	CG	ARG A 243	-61.495	-3.134	1.747	1.00	42.18	C
ATOM	3534	CD	ARG A 243	-62.826	-2.623	1.169	1.00	45.18	C
ATOM	3537	NE	ARG A 243	-63.506	-3.656	.369	1.00	48.22	N
ATOM	3539	CZ	ARG A 243	-64.374	-3.431	-.629	1.00	49.95	C
ATOM	3540	NH1	ARG A 243	-64.717	-2.182	-.991	1.00	50.13	N

ATOM	3543	NH2	ARG A 243	-64.911	-4.477	-1.273	1.00	49.77	N
ATOM	3546	C	ARG A 243	-58.438	-4.852	.068	1.00	40.08	C
ATOM	3547	O	ARG A 243	-58.084	-5.870	.665	1.00	40.45	O
ATOM	3549	N	ARG A 244	-58.236	-4.663	-1.236	1.00	40.42	N
ATOM	3550	CA	ARG A 244	-57.499	-5.621	-2.063	1.00	40.65	C
ATOM	3552	CB	ARG A 244	-57.370	-5.100	-3.503	1.00	41.29	C
ATOM	3555	CG	ARG A 244	-56.939	-6.126	-4.560	1.00	43.27	C
ATOM	3558	CD	ARG A 244	-58.090	-7.062	-4.952	1.00	46.01	C
ATOM	3561	NE	ARG A 244	-57.595	-8.337	-5.485	1.00	48.77	N
ATOM	3563	CZ	ARG A 244	-57.075	-9.332	-4.752	1.00	51.17	C
ATOM	3564	NH1	ARG A 244	-56.968	-9.234	-3.422	1.00	51.63	N
ATOM	3567	NH2	ARG A 244	-56.656	-10.448	-5.353	1.00	52.11	N
ATOM	3570	C	ARG A 244	-56.120	-5.861	-1.465	1.00	40.12	C
ATOM	3571	O	ARG A 244	-55.728	-7.000	-1.235	1.00	39.96	O
ATOM	3573	N	VAL A 245	-55.399	-4.778	-1.201	1.00	39.78	N
ATOM	3574	CA	VAL A 245	-54.099	-4.858	-.543	1.00	39.68	C
ATOM	3576	CB	VAL A 245	-53.430	-3.458	-.437	1.00	39.70	C
ATOM	3578	CG1	VAL A 245	-53.133	-2.914	-1.824	1.00	39.54	C
ATOM	3582	CG2	VAL A 245	-52.145	-3.507	.395	1.00	39.41	C
ATOM	3586	C	VAL A 245	-54.262	-5.472	.841	1.00	39.73	C
ATOM	3587	O	VAL A 245	-53.455	-6.291	1.253	1.00	39.53	O
ATOM	3589	N	GLY A 246	-55.312	-5.051	1.544	1.00	40.11	N
ATOM	3590	CA	GLY A 246	-55.705	-5.617	2.839	1.00	40.35	C
ATOM	3593	C	GLY A 246	-54.621	-5.648	3.894	1.00	40.58	C
ATOM	3594	O	GLY A 246	-54.396	-6.676	4.510	1.00	40.68	O
ATOM	3596	N	LEU A 247	-53.963	-4.522	4.129	1.00	41.15	N
ATOM	3597	CA	LEU A 247	-52.778	-4.519	4.982	1.00	41.71	C
ATOM	3599	CB	LEU A 247	-51.808	-3.422	4.537	1.00	41.60	C
ATOM	3602	CG	LEU A 247	-50.334	-3.815	4.441	1.00	40.83	C
ATOM	3604	CD1	LEU A 247	-50.158	-5.014	3.539	1.00	40.54	C
ATOM	3608	CD2	LEU A 247	-49.534	-2.644	3.918	1.00	39.92	C
ATOM	3612	C	LEU A 247	-53.136	-4.385	6.471	1.00	42.60	C
ATOM	3613	O	LEU A 247	-52.588	-5.111	7.313	1.00	42.06	O
ATOM	3615	N	ALA A 248	-54.064	-3.475	6.788	1.00	43.75	N
ATOM	3616	CA	ALA A 248	-54.603	-3.352	8.162	1.00	44.52	C
ATOM	3618	CB	ALA A 248	-55.712	-2.296	8.210	1.00	44.31	C
ATOM	3622	C	ALA A 248	-55.129	-4.706	8.695	1.00	45.10	C
ATOM	3623	O	ALA A 248	-54.969	-5.041	9.875	1.00	44.96	O
ATOM	3625	N	THR A 249	-55.742	-5.477	7.802	1.00	45.91	N
ATOM	3626	CA	THR A 249	-56.328	-6.768	8.154	1.00	46.50	C
ATOM	3628	CB	THR A 249	-57.303	-7.277	7.041	1.00	46.48	C
ATOM	3630	OG1	THR A 249	-56.583	-8.039	6.063	1.00	46.20	O
ATOM	3632	CG2	THR A 249	-58.033	-6.102	6.350	1.00	46.70	C
ATOM	3636	C	THR A 249	-55.257	-7.840	8.447	1.00	47.06	C
ATOM	3637	O	THR A 249	-55.531	-8.813	9.159	1.00	47.22	O
ATOM	3639	N	LYS A 250	-54.049	-7.661	7.905	1.00	47.52	N
ATOM	3640	CA	LYS A 250	-52.952	-8.620	8.096	1.00	47.91	C
ATOM	3642	CB	LYS A 250	-52.319	-8.976	6.744	1.00	48.09	C
ATOM	3645	CG	LYS A 250	-52.911	-10.229	6.088	1.00	48.91	C
ATOM	3648	CD	LYS A 250	-52.210	-11.510	6.585	1.00	49.88	C
ATOM	3651	CE	LYS A 250	-53.066	-12.760	6.349	1.00	50.12	C
ATOM	3654	NZ	LYS A 250	-52.385	-14.011	6.791	1.00	49.99	N
ATOM	3658	C	LYS A 250	-51.880	-8.129	9.078	1.00	48.13	C
ATOM	3659	O	LYS A 250	-51.252	-8.939	9.759	1.00	47.89	O
ATOM	3661	N	LEU A 251	-51.667	-6.812	9.134	1.00	48.63	N
ATOM	3662	CA	LEU A 251	-50.758	-6.187	10.114	1.00	48.83	C
ATOM	3664	CB	LEU A 251	-49.981	-5.020	9.485	1.00	48.74	C
ATOM	3667	CG	LEU A 251	-48.569	-5.296	8.960	1.00	48.51	C
ATOM	3669	CD1	LEU A 251	-48.472	-6.594	8.174	1.00	47.96	C
ATOM	3673	CD2	LEU A 251	-48.102	-4.115	8.116	1.00	48.76	C

ATOM	3677	C	LEU A 251	-51.559	-5.703	11.324	1.00	49.17	C
ATOM	3678	O	LEU A 251	-52.176	-4.632	11.304	1.00	49.28	O
ATOM	3680	N	HIS A 252	-51.521	-6.486	12.393	1.00	49.53	N
ATOM	3681	CA	HIS A 252	-52.470	-6.322	13.494	1.00	50.02	C
ATOM	3683	CB	HIS A 252	-52.598	-7.647	14.266	1.00	50.42	C
ATOM	3686	CG	HIS A 252	-52.860	-8.834	13.380	1.00	52.23	C
ATOM	3687	ND1	HIS A 252	-54.032	-8.989	12.667	1.00	53.75	N
ATOM	3689	CE1	HIS A 252	-53.978	-10.113	11.973	1.00	54.43	C
ATOM	3691	NE2	HIS A 252	-52.810	-10.689	12.201	1.00	54.52	N
ATOM	3693	CD2	HIS A 252	-52.090	-9.910	13.077	1.00	53.70	C
ATOM	3695	C	HIS A 252	-52.161	-5.147	14.441	1.00	49.64	C
ATOM	3696	O	HIS A 252	-52.951	-4.849	15.336	1.00	49.61	O
ATOM	3698	N	PHE A 253	-51.027	-4.483	14.229	1.00	49.37	N
ATOM	3699	CA	PHE A 253	-50.652	-3.273	14.973	1.00	49.15	C
ATOM	3701	CB	PHE A 253	-49.144	-3.273	15.260	1.00	49.07	C
ATOM	3704	CG	PHE A 253	-48.307	-3.112	14.025	1.00	48.30	C
ATOM	3705	CD1	PHE A 253	-48.028	-1.851	13.518	1.00	48.10	C
ATOM	3707	CE1	PHE A 253	-47.290	-1.702	12.349	1.00	48.21	C
ATOM	3709	CZ	PHE A 253	-46.825	-2.825	11.675	1.00	47.88	C
ATOM	3711	CE2	PHE A 253	-47.101	-4.090	12.175	1.00	47.66	C
ATOM	3713	CD2	PHE A 253	-47.842	-4.226	13.338	1.00	47.74	C
ATOM	3715	C	PHE A 253	-50.974	-2.006	14.176	1.00	49.31	C
ATOM	3716	O	PHE A 253	-50.846	-.893	14.702	1.00	49.01	O
ATOM	3718	N	ALA A 254	-51.365	-2.185	12.910	1.00	49.48	N
ATOM	3719	CA	ALA A 254	-51.412	-1.095	11.931	1.00	49.63	C
ATOM	3721	CB	ALA A 254	-51.344	-1.663	10.524	1.00	49.64	C
ATOM	3725	C	ALA A 254	-52.650	-.224	12.061	1.00	49.77	C
ATOM	3726	O	ALA A 254	-53.761	-.735	12.219	1.00	49.86	O
ATOM	3728	N	ARG A 255	-52.452	1.090	11.975	1.00	49.93	N
ATOM	3729	CA	ARG A 255	-53.562	2.042	11.908	1.00	50.27	C
ATOM	3731	CB	ARG A 255	-53.094	3.468	12.237	1.00	50.43	C
ATOM	3734	CG	ARG A 255	-52.678	3.736	13.696	1.00	50.82	C
ATOM	3737	CD	ARG A 255	-52.242	5.211	13.887	1.00	51.32	C
ATOM	3740	NE	ARG A 255	-51.003	5.522	13.155	1.00	51.81	N
ATOM	3742	CZ	ARG A 255	-50.571	6.746	12.829	1.00	51.67	C
ATOM	3743	NH1	ARG A 255	-51.264	7.839	13.146	1.00	51.63	N
ATOM	3746	NH2	ARG A 255	-49.427	6.877	12.162	1.00	51.52	N
ATOM	3749	C	ARG A 255	-54.160	2.040	10.497	1.00	50.32	C
ATOM	3750	O	ARG A 255	-53.416	2.081	9.504	1.00	50.35	O
ATOM	3752	N	ASP A 256	-55.492	2.005	10.412	1.00	50.17	N
ATOM	3753	CA	ASP A 256	-56.190	2.090	9.128	1.00	50.18	C
ATOM	3755	CB	ASP A 256	-57.259	1.009	9.037	1.00	50.37	C
ATOM	3758	CG	ASP A 256	-58.095	1.130	7.780	1.00	50.96	C
ATOM	3759	OD1	ASP A 256	-59.252	1.601	7.881	1.00	52.00	O
ATOM	3760	OD2	ASP A 256	-57.581	.787	6.693	1.00	50.94	O
ATOM	3761	C	ASP A 256	-56.838	3.470	8.945	1.00	49.85	C
ATOM	3762	O	ASP A 256	-57.812	3.804	9.624	1.00	50.13	O
ATOM	3764	N	ARG A 257	-56.315	4.263	8.013	1.00	49.11	N
ATOM	3765	CA	ARG A 257	-56.712	5.665	7.913	1.00	48.59	C
ATOM	3767	CB	ARG A 257	-55.622	6.553	8.529	1.00	48.79	C
ATOM	3770	CG	ARG A 257	-55.318	6.275	10.006	1.00	49.85	C
ATOM	3773	CD	ARG A 257	-56.538	6.480	10.908	1.00	51.01	C
ATOM	3776	NE	ARG A 257	-56.198	7.274	12.089	1.00	51.80	N
ATOM	3778	CZ	ARG A 257	-55.779	6.791	13.260	1.00	52.00	C
ATOM	3779	NH1	ARG A 257	-55.643	5.486	13.468	1.00	51.64	N
ATOM	3782	NH2	ARG A 257	-55.501	7.638	14.245	1.00	52.62	N
ATOM	3785	C	ARG A 257	-56.988	6.112	6.477	1.00	47.58	C
ATOM	3786	O	ARG A 257	-56.476	7.144	6.031	1.00	47.40	O
ATOM	3788	N	LEU A 258	-57.814	5.352	5.761	1.00	46.18	N
ATOM	3789	CA	LEU A 258	-58.094	5.672	4.358	1.00	44.78	C

ATOM	3791	CB	LEU A 258	-58.777	4.508	3.629	1.00	44.57	C
ATOM	3794	CG	LEU A 258	-58.917	4.731	2.121	1.00	43.79	C
ATOM	3796	CD1	LEU A 258	-57.557	4.661	1.478	1.00	43.89	C
ATOM	3800	CD2	LEU A 258	-59.850	3.733	1.492	1.00	42.88	C
ATOM	3804	C	LEU A 258	-58.970	6.907	4.259	1.00	43.52	C
ATOM	3805	O	LEU A 258	-58.643	7.849	3.545	1.00	43.42	O
ATOM	3807	N	ILE A 259	-60.080	6.895	4.983	1.00	41.98	N
ATOM	3808	CA	ILE A 259	-61.072	7.948	4.843	1.00	40.89	C
ATOM	3810	CB	ILE A 259	-62.349	7.703	5.693	1.00	41.02	C
ATOM	3812	CG1	ILE A 259	-62.846	6.248	5.570	1.00	41.82	C
ATOM	3815	CD1	ILE A 259	-64.029	5.872	6.508	1.00	42.29	C
ATOM	3819	CG2	ILE A 259	-63.446	8.662	5.262	1.00	40.47	C
ATOM	3823	C	ILE A 259	-60.424	9.253	5.273	1.00	39.74	C
ATOM	3824	O	ILE A 259	-60.607	10.291	4.635	1.00	39.25	O
ATOM	3826	N	GLU A 260	-59.650	9.191	6.355	1.00	38.46	N
ATOM	3827	CA	GLU A 260	-58.951	10.366	6.846	1.00	37.44	C
ATOM	3829	CB	GLU A 260	-58.183	10.072	8.140	1.00	37.74	C
ATOM	3832	CG	GLU A 260	-59.041	9.934	9.393	1.00	38.71	C
ATOM	3835	CD	GLU A 260	-59.469	8.501	9.693	1.00	41.00	C
ATOM	3836	OE1	GLU A 260	-59.273	7.594	8.842	1.00	42.70	O
ATOM	3837	OE2	GLU A 260	-60.012	8.281	10.800	1.00	42.43	O
ATOM	3838	C	GLU A 260	-57.995	10.841	5.764	1.00	35.94	C
ATOM	3839	O	GLU A 260	-57.981	12.020	5.418	1.00	36.10	O
ATOM	3841	N	SER A 261	-57.222	9.908	5.212	1.00	34.02	N
ATOM	3842	CA	SER A 261	-56.263	10.229	4.152	1.00	32.47	C
ATOM	3844	CB	SER A 261	-55.407	9.015	3.774	1.00	32.51	C
ATOM	3847	OG	SER A 261	-54.253	8.940	4.584	1.00	33.14	O
ATOM	3849	C	SER A 261	-56.894	10.764	2.889	1.00	30.82	C
ATOM	3850	O	SER A 261	-56.199	11.306	2.062	1.00	30.79	O
ATOM	3852	N	PHE A 262	-58.194	10.589	2.713	1.00	29.25	N
ATOM	3853	CA	PHE A 262	-58.861	11.091	1.514	1.00	27.98	C
ATOM	3855	CB	PHE A 262	-60.011	10.185	1.119	1.00	27.57	C
ATOM	3858	CG	PHE A 262	-60.473	10.423	-.251	1.00	26.11	C
ATOM	3859	CD1	PHE A 262	-59.763	9.914	-1.318	1.00	25.14	C
ATOM	3861	CE1	PHE A 262	-60.169	10.147	-2.605	1.00	24.32	C
ATOM	3863	CZ	PHE A 262	-61.284	10.917	-2.840	1.00	24.36	C
ATOM	3865	CE2	PHE A 262	-61.987	11.445	-1.785	1.00	25.23	C
ATOM	3867	CD2	PHE A 262	-61.575	11.201	-.492	1.00	25.61	C
ATOM	3869	C	PHE A 262	-59.394	12.507	1.707	1.00	27.29	C
ATOM	3870	O	PHE A 262	-59.275	13.359	.821	1.00	27.39	O
ATOM	3872	N	TYR A 263	-60.025	12.723	2.856	1.00	26.19	N
ATOM	3873	CA	TYR A 263	-60.415	14.052	3.310	1.00	25.22	C
ATOM	3875	CB	TYR A 263	-60.975	13.942	4.735	1.00	25.21	C
ATOM	3878	CG	TYR A 263	-61.037	15.199	5.578	1.00	25.37	C
ATOM	3879	CD1	TYR A 263	-62.020	16.152	5.385	1.00	24.95	C
ATOM	3881	CE1	TYR A 263	-62.084	17.282	6.197	1.00	26.54	C
ATOM	3883	CZ	TYR A 263	-61.158	17.455	7.229	1.00	27.44	C
ATOM	3884	OH	TYR A 263	-61.186	18.567	8.059	1.00	28.97	O
ATOM	3886	CE2	TYR A 263	-60.185	16.512	7.438	1.00	27.10	C
ATOM	3888	CD2	TYR A 263	-60.136	15.390	6.626	1.00	26.91	C
ATOM	3890	C	TYR A 263	-59.190	14.946	3.238	1.00	24.29	C
ATOM	3891	O	TYR A 263	-59.267	16.077	2.757	1.00	24.13	O
ATOM	3893	N	TRP A 264	-58.055	14.415	3.682	1.00	23.01	N
ATOM	3894	CA	TRP A 264	-56.789	15.121	3.567	1.00	22.32	C
ATOM	3896	CB	TRP A 264	-55.642	14.261	4.119	1.00	22.26	C
ATOM	3899	CG	TRP A 264	-54.326	14.874	3.860	1.00	22.16	C
ATOM	3900	CD1	TRP A 264	-53.614	14.821	2.699	1.00	22.68	C
ATOM	3902	NE1	TRP A 264	-52.456	15.538	2.818	1.00	23.08	N
ATOM	3904	CE2	TRP A 264	-52.407	16.083	4.072	1.00	23.07	C
ATOM	3905	CD2	TRP A 264	-53.579	15.686	4.753	1.00	22.18	C

ATOM	3906	CE3 TRP A 264	-53.782	16.119	6.066	1.00	22.57	C
ATOM	3908	CZ3 TRP A 264	-52.824	16.925	6.657	1.00	23.43	C
ATOM	3910	CH2 TRP A 264	-51.656	17.297	5.955	1.00	24.14	C
ATOM	3912	CZ2 TRP A 264	-51.433	16.886	4.663	1.00	23.64	C
ATOM	3914	C TRP A 264	-56.485	15.523	2.112	1.00	21.50	C
ATOM	3915	O TRP A 264	-56.179	16.682	1.820	1.00	21.29	O
ATOM	3917	N ALA A 265	-56.561	14.547	1.213	1.00	20.41	N
ATOM	3918	CA ALA A 265	-56.242	14.760	-.187	1.00	19.63	C
ATOM	3920	CB ALA A 265	-56.298	13.437	-.953	1.00	19.39	C
ATOM	3924	C ALA A 265	-57.176	15.792	-.798	1.00	19.07	C
ATOM	3925	O ALA A 265	-56.760	16.573	-1.633	1.00	18.99	O
ATOM	3927	N VAL A 266	-58.431	15.814	-.359	1.00	18.78	N
ATOM	3928	CA VAL A 266	-59.396	16.828	-.818	1.00	18.56	C
ATOM	3930	CB VAL A 266	-60.839	16.541	-.278	1.00	18.41	C
ATOM	3932	CG1 VAL A 266	-61.405	15.298	-.935	1.00	18.88	C
ATOM	3936	CG2 VAL A 266	-61.780	17.709	-.505	1.00	17.21	C
ATOM	3940	C VAL A 266	-58.947	18.257	-.450	1.00	18.43	C
ATOM	3941	O VAL A 266	-59.250	19.200	-1.171	1.00	18.83	O
ATOM	3943	N GLY A 267	-58.230	18.411	.663	1.00	17.96	N
ATOM	3944	CA GLY A 267	-57.714	19.712	1.067	1.00	17.45	C
ATOM	3947	C GLY A 267	-56.687	20.180	.070	1.00	17.20	C
ATOM	3948	O GLY A 267	-56.574	21.365	-.215	1.00	17.41	O
ATOM	3950	N VAL A 268	-55.959	19.219	-.480	1.00	16.88	N
ATOM	3951	CA VAL A 268	-54.822	19.487	-1.319	1.00	16.62	C
ATOM	3953	CB VAL A 268	-53.833	18.300	-1.290	1.00	16.19	C
ATOM	3955	CG1 VAL A 268	-52.691	18.536	-2.240	1.00	15.50	C
ATOM	3959	CG2 VAL A 268	-53.310	18.121	.103	1.00	15.30	C
ATOM	3963	C VAL A 268	-55.262	19.788	-2.736	1.00	17.15	C
ATOM	3964	O VAL A 268	-54.748	20.710	-3.350	1.00	17.12	O
ATOM	3966	N ALA A 269	-56.217	19.014	-3.243	1.00	17.99	N
ATOM	3967	CA ALA A 269	-56.688	19.158	-4.622	1.00	18.79	C
ATOM	3969	CB ALA A 269	-55.958	18.183	-5.551	1.00	18.26	C
ATOM	3973	C ALA A 269	-58.196	18.938	-4.672	1.00	19.70	C
ATOM	3974	O ALA A 269	-58.665	17.834	-4.958	1.00	19.50	O
ATOM	3976	N PHE A 270	-58.945	20.011	-4.417	1.00	21.09	N
ATOM	3977	CA PHE A 270	-60.393	19.910	-4.211	1.00	22.32	C
ATOM	3979	CB PHE A 270	-60.925	21.108	-3.404	1.00	22.49	C
ATOM	3982	CG PHE A 270	-61.193	22.321	-4.246	1.00	23.31	C
ATOM	3983	CD1 PHE A 270	-62.423	22.483	-4.883	1.00	24.16	C
ATOM	3985	CE1 PHE A 270	-62.662	23.569	-5.683	1.00	23.84	C
ATOM	3987	CZ PHE A 270	-61.675	24.497	-5.871	1.00	23.93	C
ATOM	3989	CE2 PHE A 270	-60.448	24.349	-5.242	1.00	23.33	C
ATOM	3991	CD2 PHE A 270	-60.213	23.268	-4.443	1.00	23.11	C
ATOM	3993	C PHE A 270	-61.173	19.820	-5.515	1.00	23.21	C
ATOM	3994	O PHE A 270	-62.142	19.069	-5.599	1.00	23.48	O
ATOM	3996	N GLU A 271	-60.771	20.599	-6.522	1.00	24.16	N
ATOM	3997	CA GLU A 271	-61.595	20.744	-7.741	1.00	25.01	C
ATOM	3999	CB GLU A 271	-61.065	21.863	-8.655	1.00	25.19	C
ATOM	4002	CG GLU A 271	-59.563	21.909	-8.794	1.00	27.45	C
ATOM	4005	CD GLU A 271	-58.857	22.910	-7.857	1.00	30.03	C
ATOM	4006	OE1 GLU A 271	-58.905	24.145	-8.147	1.00	29.41	O
ATOM	4007	OE2 GLU A 271	-58.227	22.431	-6.865	1.00	30.71	O
ATOM	4008	C GLU A 271	-61.799	19.391	-8.484	1.00	24.95	C
ATOM	4009	O GLU A 271	-60.972	18.480	-8.352	1.00	25.70	O
ATOM	4011	N PRO A 272	-62.918	19.242	-9.224	1.00	24.59	N
ATOM	4012	CA PRO A 272	-63.399	17.915	-9.656	1.00	24.38	C
ATOM	4014	CB PRO A 272	-64.675	18.243	-10.430	1.00	24.25	C
ATOM	4017	CG PRO A 272	-65.105	19.536	-9.869	1.00	24.53	C
ATOM	4020	CD PRO A 272	-63.855	20.295	-9.639	1.00	24.42	C
ATOM	4023	C PRO A 272	-62.475	17.111	-10.549	1.00	24.29	C

ATOM	4024	O	PRO A 272	-62.406	15.886	-10.410	1.00	24.45	O
ATOM	4025	N	GLN A 273	-61.777	17.781	-11.462	1.00	24.13	N
ATOM	4026	CA	GLN A 273	-60.998	17.064	-12.461	1.00	24.13	C
ATOM	4028	CB	GLN A 273	-60.523	17.983	-13.567	1.00	23.84	C
ATOM	4031	CG	GLN A 273	-59.556	19.041	-13.117	1.00	24.37	C
ATOM	4034	CD	GLN A 273	-60.222	20.360	-12.775	1.00	25.22	C
ATOM	4035	OE1	GLN A 273	-61.392	20.412	-12.367	1.00	25.71	O
ATOM	4036	NE2	GLN A 273	-59.468	21.441	-12.932	1.00	24.95	N
ATOM	4039	C	GLN A 273	-59.808	16.327	-11.877	1.00	24.33	C
ATOM	4040	O	GLN A 273	-59.224	15.486	-12.555	1.00	24.91	O
ATOM	4042	N	TYR A 274	-59.462	16.607	-10.626	1.00	24.27	N
ATOM	4043	CA	TYR A 274	-58.211	16.119	-10.067	1.00	24.38	C
ATOM	4045	CB	TYR A 274	-57.639	17.180	-9.138	1.00	24.36	C
ATOM	4048	CG	TYR A 274	-57.066	18.398	-9.819	1.00	24.35	C
ATOM	4049	CD1	TYR A 274	-56.249	18.291	-10.935	1.00	24.13	C
ATOM	4051	CE1	TYR A 274	-55.716	19.409	-11.542	1.00	23.75	C
ATOM	4053	CZ	TYR A 274	-55.968	20.646	-11.019	1.00	24.14	C
ATOM	4054	OH	TYR A 274	-55.431	21.762	-11.597	1.00	24.45	O
ATOM	4056	CE2	TYR A 274	-56.755	20.780	-9.903	1.00	24.92	C
ATOM	4058	CD2	TYR A 274	-57.297	19.658	-9.308	1.00	24.74	C
ATOM	4060	C	TYR A 274	-58.319	14.796	-9.310	1.00	24.63	C
ATOM	4061	O	TYR A 274	-57.681	14.613	-8.276	1.00	24.47	O
ATOM	4063	N	SER A 275	-59.097	13.852	-9.821	1.00	25.00	N
ATOM	4064	CA	SER A 275	-59.254	12.575	-9.116	1.00	25.11	C
ATOM	4066	CB	SER A 275	-60.244	11.664	-9.842	1.00	25.14	C
ATOM	4069	OG	SER A 275	-61.537	12.252	-9.852	1.00	25.96	O
ATOM	4071	C	SER A 275	-57.907	11.885	-8.963	1.00	24.94	C
ATOM	4072	O	SER A 275	-57.544	11.439	-7.877	1.00	24.57	O
ATOM	4074	N	ASP A 276	-57.149	11.822	-10.052	1.00	25.01	N
ATOM	4075	CA	ASP A 276	-55.840	11.187	-9.990	1.00	24.92	C
ATOM	4077	CB	ASP A 276	-55.102	11.257	-11.324	1.00	25.02	C
ATOM	4080	CG	ASP A 276	-55.827	10.495	-12.407	1.00	26.27	C
ATOM	4081	OD1	ASP A 276	-56.372	9.411	-12.080	1.00	26.65	O
ATOM	4082	OD2	ASP A 276	-55.880	10.991	-13.563	1.00	28.21	O
ATOM	4083	C	ASP A 276	-55.026	11.798	-8.881	1.00	24.30	C
ATOM	4084	O	ASP A 276	-54.454	11.064	-8.099	1.00	24.64	O
ATOM	4086	N	CYS A 277	-55.005	13.124	-8.772	1.00	23.61	N
ATOM	4087	CA	CYS A 277	-54.210	13.746	-7.715	1.00	22.91	C
ATOM	4089	CB	CYS A 277	-54.262	15.268	-7.748	1.00	22.83	C
ATOM	4092	SG	CYS A 277	-53.048	15.994	-6.615	1.00	22.43	S
ATOM	4094	C	CYS A 277	-54.653	13.251	-6.357	1.00	22.39	C
ATOM	4095	O	CYS A 277	-53.822	12.833	-5.559	1.00	22.23	O
ATOM	4097	N	ARG A 278	-55.962	13.262	-6.115	1.00	22.01	N
ATOM	4098	CA	ARG A 278	-56.519	12.813	-4.829	1.00	21.40	C
ATOM	4100	CB	ARG A 278	-58.021	13.053	-4.753	1.00	21.07	C
ATOM	4103	CG	ARG A 278	-58.382	14.509	-4.754	1.00	20.39	C
ATOM	4106	CD	ARG A 278	-59.852	14.698	-4.527	1.00	20.05	C
ATOM	4109	NE	ARG A 278	-60.683	14.331	-5.675	1.00	18.86	N
ATOM	4111	CZ	ARG A 278	-60.993	15.144	-6.679	1.00	18.58	C
ATOM	4112	NH1	ARG A 278	-60.530	16.391	-6.730	1.00	18.45	N
ATOM	4115	NH2	ARG A 278	-61.775	14.702	-7.647	1.00	19.04	N
ATOM	4118	C	ARG A 278	-56.224	11.358	-4.530	1.00	21.12	C
ATOM	4119	O	ARG A 278	-55.804	11.050	-3.434	1.00	21.59	O
ATOM	4121	N	ASN A 279	-56.412	10.467	-5.493	1.00	20.94	N
ATOM	4122	CA	ASN A 279	-56.168	9.044	-5.243	1.00	21.08	C
ATOM	4124	CB	ASN A 279	-56.672	8.182	-6.401	1.00	21.54	C
ATOM	4127	CG	ASN A 279	-58.199	8.346	-6.642	1.00	23.94	C
ATOM	4128	OD1	ASN A 279	-58.920	8.941	-5.823	1.00	26.47	O
ATOM	4129	ND2	ASN A 279	-58.685	7.825	-7.769	1.00	25.52	N
ATOM	4132	C	ASN A 279	-54.699	8.798	-4.972	1.00	20.30	C

ATOM	4133	O	ASN A 279	-54.337	8.200	-3.977	1.00	19.95	O
ATOM	4135	N	SER A 280	-53.851	9.297	-5.856	1.00	20.02	N
ATOM	4136	CA	SER A 280	-52.409	9.335	-5.612	1.00	19.32	C
ATOM	4138	CB	SER A 280	-51.725	10.304	-6.580	1.00	19.35	C
ATOM	4141	OG	SER A 280	-50.624	9.705	-7.212	1.00	19.77	O
ATOM	4143	C	SER A 280	-52.122	9.745	-4.166	1.00	18.55	C
ATOM	4144	O	SER A 280	-51.514	8.976	-3.436	1.00	18.98	O
ATOM	4146	N	VAL A 281	-52.580	10.923	-3.743	1.00	17.45	N
ATOM	4147	CA	VAL A 281	-52.235	11.444	-2.416	1.00	17.01	C
ATOM	4149	CB	VAL A 281	-52.650	12.934	-2.235	1.00	17.20	C
ATOM	4151	CG1	VAL A 281	-52.492	13.376	-.787	1.00	16.76	C
ATOM	4155	CG2	VAL A 281	-51.826	13.847	-3.128	1.00	17.13	C
ATOM	4159	C	VAL A 281	-52.861	10.607	-1.301	1.00	16.82	C
ATOM	4160	O	VAL A 281	-52.217	10.314	-.289	1.00	16.81	O
ATOM	4162	N	ALA A 282	-54.118	10.225	-1.481	1.00	16.55	N
ATOM	4163	CA	ALA A 282	-54.774	9.304	-.559	1.00	16.31	C
ATOM	4165	CB	ALA A 282	-56.174	8.951	-1.060	1.00	16.08	C
ATOM	4169	C	ALA A 282	-53.946	8.036	-.355	1.00	16.31	C
ATOM	4170	O	ALA A 282	-53.671	7.653	.762	1.00	16.16	O
ATOM	4172	N	LYS A 283	-53.539	7.395	-1.443	1.00	16.92	N
ATOM	4173	CA	LYS A 283	-52.773	6.139	-1.368	1.00	17.39	C
ATOM	4175	CB	LYS A 283	-52.454	5.596	-2.773	1.00	17.31	C
ATOM	4178	CG	LYS A 283	-53.680	5.066	-3.519	1.00	17.85	C
ATOM	4181	CD	LYS A 283	-53.389	4.527	-4.941	1.00	18.98	C
ATOM	4184	CE	LYS A 283	-54.709	4.275	-5.718	1.00	19.79	C
ATOM	4187	NZ	LYS A 283	-54.695	3.091	-6.625	1.00	20.11	N
ATOM	4191	C	LYS A 283	-51.483	6.338	-.591	1.00	17.82	C
ATOM	4192	O	LYS A 283	-51.170	5.577	.331	1.00	17.37	O
ATOM	4194	N	MET A 284	-50.744	7.381	-.964	1.00	18.49	N
ATOM	4195	CA	MET A 284	-49.434	7.622	-.373	1.00	18.82	C
ATOM	4197	CB	MET A 284	-48.688	8.795	-1.047	1.00	18.84	C
ATOM	4200	CG	MET A 284	-48.118	8.506	-2.444	1.00	18.47	C
ATOM	4203	SD	MET A 284	-47.540	6.819	-2.718	1.00	19.50	S
ATOM	4204	CE	MET A 284	-49.088	5.982	-3.094	1.00	19.36	C
ATOM	4208	C	MET A 284	-49.575	7.848	1.115	1.00	18.95	C
ATOM	4209	O	MET A 284	-48.901	7.166	1.892	1.00	19.29	O
ATOM	4211	N	PHE A 285	-50.465	8.762	1.507	1.00	18.94	N
ATOM	4212	CA	PHE A 285	-50.671	9.082	2.930	1.00	19.26	C
ATOM	4214	CB	PHE A 285	-51.714	10.185	3.073	1.00	19.72	C
ATOM	4217	CG	PHE A 285	-51.575	11.016	4.328	1.00	21.97	C
ATOM	4218	CD1	PHE A 285	-50.335	11.202	4.954	1.00	24.49	C
ATOM	4220	CE1	PHE A 285	-50.209	11.996	6.086	1.00	25.11	C
ATOM	4222	CZ	PHE A 285	-51.321	12.634	6.597	1.00	26.40	C
ATOM	4224	CE2	PHE A 285	-52.566	12.470	5.977	1.00	26.44	C
ATOM	4226	CD2	PHE A 285	-52.678	11.670	4.846	1.00	24.55	C
ATOM	4228	C	PHE A 285	-51.097	7.879	3.770	1.00	18.74	C
ATOM	4229	O	PHE A 285	-50.700	7.729	4.924	1.00	18.31	O
ATOM	4231	N	SER A 286	-51.903	7.019	3.169	1.00	18.49	N
ATOM	4232	CA	SER A 286	-52.281	5.772	3.791	1.00	18.38	C
ATOM	4234	CB	SER A 286	-53.302	5.065	2.926	1.00	18.45	C
ATOM	4237	OG	SER A 286	-54.389	5.949	2.687	1.00	19.82	O
ATOM	4239	C	SER A 286	-51.070	4.888	4.048	1.00	18.08	C
ATOM	4240	O	SER A 286	-50.931	4.376	5.146	1.00	18.18	O
ATOM	4242	N	PHE A 287	-50.195	4.714	3.054	1.00	17.72	N
ATOM	4243	CA	PHE A 287	-48.911	4.024	3.279	1.00	17.23	C
ATOM	4245	CB	PHE A 287	-48.133	3.801	1.977	1.00	17.09	C
ATOM	4248	CG	PHE A 287	-48.513	2.555	1.258	1.00	17.26	C
ATOM	4249	CD1	PHE A 287	-48.220	1.315	1.802	1.00	18.31	C
ATOM	4251	CE1	PHE A 287	-48.583	.129	1.131	1.00	18.47	C
ATOM	4253	CZ	PHE A 287	-49.237	.198	-.091	1.00	17.88	C

ATOM	4255	CE2 PHE A 287	-49.530	1.441	-.642	1.00	17.57	C
ATOM	4257	CD2 PHE A 287	-49.168	2.608	.031	1.00	17.79	C
ATOM	4259	C PHE A 287	-48.022	4.787	4.277	1.00	16.92	C
ATOM	4260	O PHE A 287	-47.369	4.173	5.143	1.00	17.18	O
ATOM	4262	N VAL A 288	-47.989	6.113	4.174	1.00	16.00	N
ATOM	4263	CA VAL A 288	-47.148	6.883	5.080	1.00	15.45	C
ATOM	4265	CB VAL A 288	-47.289	8.399	4.863	1.00	15.21	C
ATOM	4267	CG1 VAL A 288	-46.433	9.162	5.852	1.00	14.24	C
ATOM	4271	CG2 VAL A 288	-46.898	8.757	3.450	1.00	14.46	C
ATOM	4275	C VAL A 288	-47.479	6.499	6.527	1.00	15.46	C
ATOM	4276	O VAL A 288	-46.590	6.160	7.299	1.00	14.79	O
ATOM	4278	N THR A 289	-48.759	6.494	6.880	1.00	15.62	N
ATOM	4279	CA THR A 289	-49.110	6.272	8.274	1.00	16.09	C
ATOM	4281	CB THR A 289	-50.602	6.557	8.588	1.00	15.72	C
ATOM	4283	OG1 THR A 289	-51.422	5.827	7.708	1.00	16.52	O
ATOM	4285	CG2 THR A 289	-50.927	8.007	8.371	1.00	16.56	C
ATOM	4289	C THR A 289	-48.677	4.879	8.727	1.00	16.34	C
ATOM	4290	O THR A 289	-48.234	4.707	9.881	1.00	16.68	O
ATOM	4292	N ILE A 290	-48.737	3.904	7.820	1.00	16.29	N
ATOM	4293	CA ILE A 290	-48.347	2.538	8.166	1.00	16.56	C
ATOM	4295	CB ILE A 290	-48.734	1.516	7.096	1.00	16.48	C
ATOM	4297	CG1 ILE A 290	-50.247	1.440	6.944	1.00	16.91	C
ATOM	4300	CD1 ILE A 290	-50.677	.564	5.811	1.00	16.52	C
ATOM	4304	CG2 ILE A 290	-48.262	.145	7.488	1.00	16.52	C
ATOM	4308	C ILE A 290	-46.842	2.456	8.402	1.00	16.94	C
ATOM	4309	O ILE A 290	-46.404	1.968	9.443	1.00	17.06	O
ATOM	4311	N ILE A 291	-46.048	2.942	7.451	1.00	17.22	N
ATOM	4312	CA ILE A 291	-44.595	2.969	7.640	1.00	17.20	C
ATOM	4314	CB ILE A 291	-43.842	3.534	6.405	1.00	17.32	C
ATOM	4316	CG1 ILE A 291	-44.172	2.744	5.125	1.00	17.61	C
ATOM	4319	CD1 ILE A 291	-44.004	1.265	5.250	1.00	18.21	C
ATOM	4323	CG2 ILE A 291	-42.336	3.550	6.637	1.00	16.64	C
ATOM	4327	C ILE A 291	-44.262	3.794	8.886	1.00	17.36	C
ATOM	4328	O ILE A 291	-43.420	3.413	9.670	1.00	17.06	O
ATOM	4330	N ASP A 292	-44.950	4.906	9.093	1.00	18.12	N
ATOM	4331	CA ASP A 292	-44.627	5.755	10.227	1.00	18.83	C
ATOM	4333	CB ASP A 292	-45.511	7.017	10.285	1.00	19.22	C
ATOM	4336	CG ASP A 292	-45.185	7.920	11.485	1.00	19.93	C
ATOM	4337	OD1 ASP A 292	-44.062	8.460	11.544	1.00	21.80	O
ATOM	4338	OD2 ASP A 292	-46.048	8.080	12.375	1.00	20.46	O
ATOM	4339	C ASP A 292	-44.745	4.947	11.509	1.00	18.96	C
ATOM	4340	O ASP A 292	-43.902	5.094	12.394	1.00	19.05	O
ATOM	4342	N ASP A 293	-45.774	4.097	11.610	1.00	19.05	N
ATOM	4343	CA ASP A 293	-45.951	3.273	12.817	1.00	18.99	C
ATOM	4345	CB ASP A 293	-47.237	2.443	12.790	1.00	18.95	C
ATOM	4348	CG ASP A 293	-48.496	3.272	12.921	1.00	19.27	C
ATOM	4349	OD1 ASP A 293	-48.483	4.361	13.523	1.00	21.51	O
ATOM	4350	OD2 ASP A 293	-49.536	2.812	12.421	1.00	19.61	O
ATOM	4351	C ASP A 293	-44.783	2.322	12.935	1.00	19.07	C
ATOM	4352	O ASP A 293	-44.244	2.135	14.024	1.00	19.35	O
ATOM	4354	N ILE A 294	-44.390	1.730	11.805	1.00	19.10	N
ATOM	4355	CA ILE A 294	-43.320	.727	11.786	1.00	19.04	C
ATOM	4357	CB ILE A 294	-43.137	.109	10.386	1.00	18.49	C
ATOM	4359	CG1 ILE A 294	-44.343	-.754	10.038	1.00	17.60	C
ATOM	4362	CD1 ILE A 294	-44.204	-1.504	8.752	1.00	16.77	C
ATOM	4366	CG2 ILE A 294	-41.895	-.745	10.348	1.00	18.77	C
ATOM	4370	C ILE A 294	-41.974	1.273	12.298	1.00	19.55	C
ATOM	4371	O ILE A 294	-41.234	.550	12.974	1.00	19.75	O
ATOM	4373	N TYR A 295	-41.649	2.525	11.983	1.00	19.75	N
ATOM	4374	CA TYR A 295	-40.411	3.103	12.476	1.00	19.91	C

ATOM	4376	CB	TYR	A	295	-39.856	4.178	11.542	1.00	19.66	C
ATOM	4379	CG	TYR	A	295	-39.206	3.695	10.245	1.00	18.19	C
ATOM	4380	CD1	TYR	A	295	-37.821	3.727	10.075	1.00	16.19	C
ATOM	4382	CE1	TYR	A	295	-37.222	3.331	8.882	1.00	14.68	C
ATOM	4384	CZ	TYR	A	295	-38.004	2.906	7.828	1.00	15.04	C
ATOM	4385	OH	TYR	A	295	-37.428	2.485	6.631	1.00	12.25	O
ATOM	4387	CE2	TYR	A	295	-39.383	2.871	7.975	1.00	16.48	C
ATOM	4389	CD2	TYR	A	295	-39.975	3.271	9.172	1.00	17.01	C
ATOM	4391	C	TYR	A	295	-40.647	3.677	13.858	1.00	20.73	C
ATOM	4392	O	TYR	A	295	-39.749	3.635	14.717	1.00	21.64	O
ATOM	4394	N	ASP	A	296	-41.843	4.200	14.103	1.00	21.38	N
ATOM	4395	CA	ASP	A	296	-42.100	4.865	15.391	1.00	22.20	C
ATOM	4397	CB	ASP	A	296	-43.421	5.644	15.391	1.00	22.69	C
ATOM	4400	CG	ASP	A	296	-43.567	6.548	16.607	1.00	23.85	C
ATOM	4401	OD1	ASP	A	296	-42.568	7.179	17.014	1.00	25.49	O
ATOM	4402	OD2	ASP	A	296	-44.689	6.648	17.145	1.00	26.12	O
ATOM	4403	C	ASP	A	296	-42.106	3.890	16.549	1.00	22.10	C
ATOM	4404	O	ASP	A	296	-41.376	4.082	17.506	1.00	22.14	O
ATOM	4406	N	VAL	A	297	-42.900	2.828	16.435	1.00	22.25	N
ATOM	4407	CA	VAL	A	297	-43.163	1.937	17.570	1.00	22.22	C
ATOM	4409	CB	VAL	A	297	-44.635	2.116	18.057	1.00	22.07	C
ATOM	4411	CG1	VAL	A	297	-44.958	3.587	18.215	1.00	21.63	C
ATOM	4415	CG2	VAL	A	297	-45.618	1.459	17.099	1.00	20.76	C
ATOM	4419	C	VAL	A	297	-42.861	.422	17.373	1.00	22.49	C
ATOM	4420	O	VAL	A	297	-42.517	-.255	18.343	1.00	22.40	O
ATOM	4422	N	TYR	A	298	-42.990	-.117	16.157	1.00	22.69	N
ATOM	4423	CA	TYR	A	298	-43.065	-1.579	15.993	1.00	22.90	C
ATOM	4425	CB	TYR	A	298	-44.089	-1.974	14.934	1.00	22.77	C
ATOM	4428	CG	TYR	A	298	-44.341	-3.469	14.932	1.00	23.74	C
ATOM	4429	CD1	TYR	A	298	-45.249	-4.043	15.819	1.00	25.72	C
ATOM	4431	CE1	TYR	A	298	-45.484	-5.435	15.833	1.00	26.12	C
ATOM	4433	CZ	TYR	A	298	-44.798	-6.247	14.954	1.00	25.59	C
ATOM	4434	OH	TYR	A	298	-45.027	-7.593	14.970	1.00	25.33	O
ATOM	4436	CE2	TYR	A	298	-43.890	-5.702	14.062	1.00	24.79	C
ATOM	4438	CD2	TYR	A	298	-43.661	-4.317	14.061	1.00	24.20	C
ATOM	4440	C	TYR	A	298	-41.750	-2.258	15.649	1.00	23.02	C
ATOM	4441	O	TYR	A	298	-41.387	-3.257	16.256	1.00	22.80	O
ATOM	4443	N	GLY	A	299	-41.069	-1.747	14.634	1.00	23.47	N
ATOM	4444	CA	GLY	A	299	-39.839	-2.361	14.149	1.00	23.36	C
ATOM	4447	C	GLY	A	299	-38.631	-1.978	14.984	1.00	23.27	C
ATOM	4448	O	GLY	A	299	-38.532	-.858	15.501	1.00	23.35	O
ATOM	4450	N	THR	A	300	-37.702	-2.918	15.098	1.00	23.11	N
ATOM	4451	CA	THR	A	300	-36.459	-2.682	15.797	1.00	22.90	C
ATOM	4453	CB	THR	A	300	-35.869	-3.957	16.386	1.00	22.83	C
ATOM	4455	OG1	THR	A	300	-35.328	-4.756	15.328	1.00	22.61	O
ATOM	4457	CG2	THR	A	300	-36.928	-4.732	17.162	1.00	22.18	C
ATOM	4461	C	THR	A	300	-35.482	-2.120	14.796	1.00	23.00	C
ATOM	4462	O	THR	A	300	-35.606	-2.364	13.602	1.00	23.13	O
ATOM	4464	N	LEU	A	301	-34.496	-1.394	15.309	1.00	23.07	N
ATOM	4465	CA	LEU	A	301	-33.596	-.572	14.498	1.00	22.84	C
ATOM	4467	CB	LEU	A	301	-32.515	.007	15.403	1.00	22.80	C
ATOM	4470	CG	LEU	A	301	-31.965	1.412	15.182	1.00	22.39	C
ATOM	4472	CD1	LEU	A	301	-32.983	2.365	14.587	1.00	22.44	C
ATOM	4476	CD2	LEU	A	301	-31.482	1.919	16.536	1.00	22.19	C
ATOM	4480	C	LEU	A	301	-32.964	-1.358	13.361	1.00	22.95	C
ATOM	4481	O	LEU	A	301	-32.910	-.888	12.232	1.00	22.50	O
ATOM	4483	N	ASP	A	302	-32.511	-2.569	13.673	1.00	23.35	N
ATOM	4484	CA	ASP	A	302	-31.929	-3.469	12.670	1.00	23.65	C
ATOM	4486	CB	ASP	A	302	-31.332	-4.747	13.328	1.00	24.12	C
ATOM	4489	CG	ASP	A	302	-29.984	-4.498	14.057	1.00	25.23	C

ATOM	4490	OD1 ASP A 302	-29.003	-4.051	13.409	1.00	26.43	O
ATOM	4491	OD2 ASP A 302	-29.896	-4.780	15.277	1.00	26.92	O
ATOM	4492	C ASP A 302	-32.945	-3.859	11.581	1.00	23.11	C
ATOM	4493	O ASP A 302	-32.544	-4.086	10.445	1.00	23.05	O
ATOM	4495	N GLU A 303	-34.235	-3.965	11.927	1.00	22.63	N
ATOM	4496	CA GLU A 303	-35.295	-4.248	10.931	1.00	22.40	C
ATOM	4498	CB GLU A 303	-36.625	-4.675	11.589	1.00	22.20	C
ATOM	4501	CG GLU A 303	-36.537	-5.956	12.426	1.00	22.66	C
ATOM	4504	CD GLU A 303	-37.832	-6.341	13.159	1.00	22.69	C
ATOM	4505	OE1 GLU A 303	-38.614	-5.449	13.568	1.00	22.23	O
ATOM	4506	OE2 GLU A 303	-38.051	-7.558	13.338	1.00	21.64	O
ATOM	4507	C GLU A 303	-35.533	-3.017	10.067	1.00	22.35	C
ATOM	4508	O GLU A 303	-35.658	-3.105	8.841	1.00	22.27	O
ATOM	4510	N LEU A 304	-35.597	-1.866	10.724	1.00	22.34	N
ATOM	4511	CA LEU A 304	-35.776	-.607	10.040	1.00	22.26	C
ATOM	4513	CB LEU A 304	-35.861	.537	11.052	1.00	22.49	C
ATOM	4516	CG LEU A 304	-37.089	.517	11.973	1.00	23.05	C
ATOM	4518	CD1 LEU A 304	-37.109	1.708	12.932	1.00	23.46	C
ATOM	4522	CD2 LEU A 304	-38.353	.511	11.142	1.00	23.96	C
ATOM	4526	C LEU A 304	-34.634	-.383	9.063	1.00	22.09	C
ATOM	4527	O LEU A 304	-34.873	.094	7.969	1.00	22.32	O
ATOM	4529	N GLU A 305	-33.408	-.740	9.450	1.00	21.91	N
ATOM	4530	CA GLU A 305	-32.259	-.685	8.541	1.00	21.93	C
ATOM	4532	CB GLU A 305	-30.988	-1.194	9.219	1.00	22.35	C
ATOM	4535	CG GLU A 305	-30.363	-.258	10.249	1.00	24.12	C
ATOM	4538	CD GLU A 305	-29.751	1.002	9.649	1.00	26.11	C
ATOM	4539	OE1 GLU A 305	-29.275	1.846	10.454	1.00	25.68	O
ATOM	4540	OE2 GLU A 305	-29.756	1.146	8.391	1.00	27.67	O
ATOM	4541	C GLU A 305	-32.475	-1.524	7.295	1.00	21.54	C
ATOM	4542	O GLU A 305	-32.010	-1.162	6.220	1.00	21.75	O
ATOM	4544	N LEU A 306	-33.155	-2.660	7.448	1.00	21.21	N
ATOM	4545	CA LEU A 306	-33.459	-3.557	6.316	1.00	20.66	C
ATOM	4547	CB LEU A 306	-33.880	-4.956	6.795	1.00	20.26	C
ATOM	4550	CG LEU A 306	-32.725	-5.854	7.220	1.00	19.57	C
ATOM	4552	CD1 LEU A 306	-33.261	-7.150	7.777	1.00	20.18	C
ATOM	4556	CD2 LEU A 306	-31.791	-6.120	6.061	1.00	18.29	C
ATOM	4560	C LEU A 306	-34.535	-2.967	5.416	1.00	20.23	C
ATOM	4561	O LEU A 306	-34.412	-2.966	4.197	1.00	19.97	O
ATOM	4563	N PHE A 307	-35.588	-2.457	6.023	1.00	19.99	N
ATOM	4564	CA PHE A 307	-36.635	-1.861	5.237	1.00	20.00	C
ATOM	4566	CB PHE A 307	-37.771	-1.388	6.118	1.00	20.08	C
ATOM	4569	CG PHE A 307	-39.011	-1.115	5.370	1.00	19.05	C
ATOM	4570	CD1 PHE A 307	-39.895	-2.118	5.119	1.00	18.77	C
ATOM	4572	CE1 PHE A 307	-41.037	-1.873	4.421	1.00	20.07	C
ATOM	4574	CZ PHE A 307	-41.296	-.621	3.959	1.00	19.63	C
ATOM	4576	CE2 PHE A 307	-40.411	.391	4.199	1.00	19.27	C
ATOM	4578	CD2 PHE A 307	-39.278	.142	4.900	1.00	19.17	C
ATOM	4580	C PHE A 307	-36.087	-.695	4.447	1.00	20.25	C
ATOM	4581	O PHE A 307	-36.356	-.589	3.242	1.00	20.08	O
ATOM	4583	N THR A 308	-35.325	.167	5.137	1.00	20.51	N
ATOM	4584	CA THR A 308	-34.737	1.379	4.538	1.00	20.60	C
ATOM	4586	CB THR A 308	-33.860	2.177	5.541	1.00	20.41	C
ATOM	4588	OG1 THR A 308	-34.672	2.661	6.611	1.00	20.05	O
ATOM	4590	CG2 THR A 308	-33.198	3.376	4.864	1.00	20.38	C
ATOM	4594	C THR A 308	-33.895	1.009	3.331	1.00	20.82	C
ATOM	4595	O THR A 308	-34.057	1.571	2.243	1.00	20.45	O
ATOM	4597	N ASP A 309	-33.020	.031	3.522	1.00	21.32	N
ATOM	4598	CA ASP A 309	-32.138	-.399	2.449	1.00	21.88	C
ATOM	4600	CB ASP A 309	-31.042	-1.324	2.971	1.00	22.44	C
ATOM	4603	CG ASP A 309	-30.164	-1.847	1.857	1.00	25.69	C

ATOM	4604	OD1 ASP A 309	-29.599	-1.010	1.089	1.00	28.02	O
ATOM	4605	OD2 ASP A 309	-30.076	-3.100	1.730	1.00	30.59	O
ATOM	4606	C ASP A 309	-32.928	-1.085	1.357	1.00	21.17	C
ATOM	4607	O ASP A 309	-32.601	-.948	.188	1.00	20.95	O
ATOM	4609	N ALA A 310	-33.965	-1.818	1.757	1.00	20.99	N
ATOM	4610	CA ALA A 310	-34.885	-2.474	.828	1.00	20.91	C
ATOM	4612	CB ALA A 310	-35.931	-3.244	1.591	1.00	20.77	C
ATOM	4616	C ALA A 310	-35.567	-1.501	-.123	1.00	21.08	C
ATOM	4617	O ALA A 310	-35.745	-1.816	-1.305	1.00	20.76	O
ATOM	4619	N VAL A 311	-35.958	-.334	.402	1.00	21.51	N
ATOM	4620	CA VAL A 311	-36.594	.719	-.396	1.00	21.50	C
ATOM	4622	CB VAL A 311	-37.424	1.715	.477	1.00	21.73	C
ATOM	4624	CG1 VAL A 311	-38.724	1.076	.941	1.00	21.01	C
ATOM	4628	CG2 VAL A 311	-37.765	2.993	-.292	1.00	21.88	C
ATOM	4632	C VAL A 311	-35.553	1.445	-1.230	1.00	21.67	C
ATOM	4633	O VAL A 311	-35.811	1.728	-2.389	1.00	21.71	O
ATOM	4635	N GLU A 312	-34.381	1.724	-.663	1.00	22.25	N
ATOM	4636	CA GLU A 312	-33.248	2.284	-1.443	1.00	23.00	C
ATOM	4638	CB GLU A 312	-31.951	2.331	-.608	1.00	23.36	C
ATOM	4641	CG GLU A 312	-31.897	3.383	.511	1.00	24.72	C
ATOM	4644	CD GLU A 312	-30.526	3.469	1.189	1.00	27.09	C
ATOM	4645	OE1 GLU A 312	-30.083	4.608	1.458	1.00	30.07	O
ATOM	4646	OE2 GLU A 312	-29.885	2.418	1.454	1.00	27.86	O
ATOM	4647	C GLU A 312	-32.954	1.498	-2.731	1.00	23.13	C
ATOM	4648	O GLU A 312	-32.851	2.064	-3.803	1.00	22.88	O
ATOM	4650	N ARG A 313	-32.819	.188	-2.615	1.00	23.74	N
ATOM	4651	CA ARG A 313	-32.400	-.628	-3.743	1.00	24.46	C
ATOM	4653	CB ARG A 313	-31.758	-1.916	-3.232	1.00	25.01	C
ATOM	4656	CG ARG A 313	-30.421	-1.667	-2.481	1.00	27.74	C
ATOM	4659	CD ARG A 313	-29.715	-2.963	-2.072	1.00	31.41	C
ATOM	4662	NE ARG A 313	-30.675	-3.921	-1.499	1.00	34.94	N
ATOM	4664	CZ ARG A 313	-31.203	-4.974	-2.138	1.00	37.87	C
ATOM	4665	NH1 ARG A 313	-30.861	-5.282	-3.402	1.00	38.35	N
ATOM	4668	NH2 ARG A 313	-32.079	-5.746	-1.495	1.00	38.82	N
ATOM	4671	C ARG A 313	-33.523	-.914	-4.739	1.00	24.27	C
ATOM	4672	O ARG A 313	-33.257	-1.178	-5.898	1.00	24.17	O
ATOM	4674	N TRP A 314	-34.770	-.873	-4.282	1.00	24.59	N
ATOM	4675	CA TRP A 314	-35.939	-1.042	-5.142	1.00	24.66	C
ATOM	4677	CB TRP A 314	-36.175	.234	-5.961	1.00	24.39	C
ATOM	4680	CG TRP A 314	-37.575	.382	-6.386	1.00	22.61	C
ATOM	4681	CD1 TRP A 314	-38.073	.204	-7.635	1.00	21.66	C
ATOM	4683	NE1 TRP A 314	-39.429	.403	-7.634	1.00	21.42	N
ATOM	4685	CE2 TRP A 314	-39.829	.707	-6.360	1.00	20.99	C
ATOM	4686	CD2 TRP A 314	-38.683	.701	-5.550	1.00	21.41	C
ATOM	4687	CE3 TRP A 314	-38.817	.991	-4.191	1.00	21.48	C
ATOM	4689	CZ3 TRP A 314	-40.080	1.277	-3.693	1.00	21.38	C
ATOM	4691	CH2 TRP A 314	-41.204	1.269	-4.523	1.00	20.84	C
ATOM	4693	CZ2 TRP A 314	-41.099	.985	-5.857	1.00	20.87	C
ATOM	4695	C TRP A 314	-35.831	-2.286	-6.038	1.00	25.65	C
ATOM	4696	O TRP A 314	-36.069	-2.238	-7.259	1.00	25.69	O
ATOM	4698	N ASP A 315	-35.480	-3.400	-5.404	1.00	26.79	N
ATOM	4699	CA ASP A 315	-35.253	-4.667	-6.087	1.00	28.01	C
ATOM	4701	CB ASP A 315	-33.801	-5.097	-5.890	1.00	28.07	C
ATOM	4704	CG ASP A 315	-33.553	-6.536	-6.296	1.00	29.92	C
ATOM	4705	OD1 ASP A 315	-34.262	-7.026	-7.197	1.00	32.58	O
ATOM	4706	OD2 ASP A 315	-32.648	-7.194	-5.722	1.00	32.73	O
ATOM	4707	C ASP A 315	-36.202	-5.710	-5.518	1.00	28.91	C
ATOM	4708	O ASP A 315	-36.019	-6.163	-4.399	1.00	29.17	O
ATOM	4710	N VAL A 316	-37.218	-6.091	-6.283	1.00	30.25	N
ATOM	4711	CA VAL A 316	-38.233	-7.012	-5.778	1.00	31.47	C

ATOM	4713	CB	VAL A 316	-39.396	-7.156	-6.755	1.00	31.56	C
ATOM	4715	CG1	VAL A 316	-40.668	-7.590	-6.027	1.00	30.99	C
ATOM	4719	CG2	VAL A 316	-39.033	-8.141	-7.862	1.00	31.84	C
ATOM	4723	C	VAL A 316	-37.663	-8.405	-5.540	1.00	32.86	C
ATOM	4724	O	VAL A 316	-38.170	-9.153	-4.708	1.00	32.99	O
ATOM	4726	N	ASN A 317	-36.607	-8.750	-6.278	1.00	34.55	N
ATOM	4727	CA	ASN A 317	-35.940	-10.060	-6.160	1.00	35.53	C
ATOM	4729	CB	ASN A 317	-35.013	-10.298	-7.367	1.00	35.68	C
ATOM	4732	CG	ASN A 317	-35.752	-10.264	-8.713	1.00	35.92	C
ATOM	4733	OD1	ASN A 317	-36.634	-11.092	-8.971	1.00	36.95	O
ATOM	4734	ND2	ASN A 317	-35.368	-9.323	-9.585	1.00	34.17	N
ATOM	4737	C	ASN A 317	-35.126	-10.225	-4.871	1.00	36.40	C
ATOM	4738	O	ASN A 317	-34.385	-11.195	-4.747	1.00	36.47	O
ATOM	4740	N	ALA A 318	-35.239	-9.265	-3.944	1.00	37.52	N
ATOM	4741	CA	ALA A 318	-34.614	-9.335	-2.613	1.00	38.45	C
ATOM	4743	CB	ALA A 318	-33.371	-8.465	-2.557	1.00	38.43	C
ATOM	4747	C	ALA A 318	-35.631	-8.902	-1.554	1.00	39.28	C
ATOM	4748	O	ALA A 318	-35.361	-8.094	-.662	1.00	39.55	O
ATOM	4750	N	ILE A 319	-36.823	-9.456	-1.698	1.00	40.17	N
ATOM	4751	CA	ILE A 319	-37.905	-9.310	-.740	1.00	40.49	C
ATOM	4753	CB	ILE A 319	-39.275	-9.588	-1.469	1.00	40.64	C
ATOM	4755	CG1	ILE A 319	-40.473	-9.078	-.683	1.00	40.88	C
ATOM	4758	CD1	ILE A 319	-41.799	-9.535	-1.280	1.00	40.86	C
ATOM	4762	CG2	ILE A 319	-39.460	-11.082	-1.801	1.00	40.50	C
ATOM	4766	C	ILE A 319	-37.656	-10.319	.393	1.00	40.69	C
ATOM	4767	O	ILE A 319	-38.136	-10.138	1.504	1.00	40.88	O
ATOM	4769	N	ASN A 320	-36.886	-11.374	.104	1.00	40.79	N
ATOM	4770	CA	ASN A 320	-36.689	-12.483	1.049	1.00	40.64	C
ATOM	4772	CB	ASN A 320	-36.240	-13.761	.314	1.00	40.73	C
ATOM	4775	CG	ASN A 320	-37.370	-14.422	-.468	1.00	40.98	C
ATOM	4776	OD1	ASN A 320	-38.556	-14.309	-.119	1.00	40.77	O
ATOM	4777	ND2	ASN A 320	-37.001	-15.131	-1.529	1.00	41.31	N
ATOM	4780	C	ASN A 320	-35.711	-12.175	2.169	1.00	40.16	C
ATOM	4781	O	ASN A 320	-35.546	-12.978	3.077	1.00	40.15	O
ATOM	4783	N	ASP A 321	-35.067	-11.017	2.101	1.00	39.66	N
ATOM	4784	CA	ASP A 321	-34.111	-10.615	3.126	1.00	39.51	C
ATOM	4786	CB	ASP A 321	-33.114	-9.571	2.575	1.00	40.17	C
ATOM	4789	CG	ASP A 321	-32.595	-9.904	1.152	1.00	42.12	C
ATOM	4790	OD1	ASP A 321	-32.425	-11.115	.820	1.00	44.15	O
ATOM	4791	OD2	ASP A 321	-32.354	-8.934	.375	1.00	43.28	O
ATOM	4792	C	ASP A 321	-34.851	-10.030	4.337	1.00	38.27	C
ATOM	4793	O	ASP A 321	-34.304	-9.987	5.443	1.00	38.15	O
ATOM	4795	N	LEU A 322	-36.089	-9.583	4.108	1.00	36.72	N
ATOM	4796	CA	LEU A 322	-36.887	-8.875	5.104	1.00	35.38	C
ATOM	4798	CB	LEU A 322	-37.894	-7.942	4.422	1.00	35.15	C
ATOM	4801	CG	LEU A 322	-37.370	-6.840	3.503	1.00	34.83	C
ATOM	4803	CD1	LEU A 322	-38.479	-6.286	2.626	1.00	34.45	C
ATOM	4807	CD2	LEU A 322	-36.742	-5.740	4.316	1.00	34.89	C
ATOM	4811	C	LEU A 322	-37.683	-9.834	5.963	1.00	34.57	C
ATOM	4812	O	LEU A 322	-37.975	-10.933	5.527	1.00	34.27	O
ATOM	4814	N	PRO A 323	-38.039	-9.402	7.189	1.00	33.93	N
ATOM	4815	CA	PRO A 323	-39.070	-9.917	8.067	1.00	33.49	C
ATOM	4817	CB	PRO A 323	-39.151	-8.840	9.141	1.00	33.31	C
ATOM	4820	CG	PRO A 323	-37.791	-8.419	9.311	1.00	33.67	C
ATOM	4823	CD	PRO A 323	-37.154	-8.496	7.941	1.00	34.20	C
ATOM	4826	C	PRO A 323	-40.429	-10.036	7.425	1.00	33.33	C
ATOM	4827	O	PRO A 323	-40.776	-9.232	6.579	1.00	33.35	O
ATOM	4828	N	ASP A 324	-41.209	-11.005	7.891	1.00	33.37	N
ATOM	4829	CA	ASP A 324	-42.511	-11.329	7.313	1.00	33.36	C
ATOM	4831	CB	ASP A 324	-43.137	-12.542	8.037	1.00	33.59	C

ATOM	4834	CG	ASP	A 324	-42.496	-13.881	7.619	1.00	34.18	C
ATOM	4835	OD1	ASP	A 324	-41.885	-13.915	6.518	1.00	36.44	O
ATOM	4836	OD2	ASP	A 324	-42.607	-14.887	8.371	1.00	32.19	O
ATOM	4837	C	ASP	A 324	-43.484	-10.149	7.289	1.00	32.90	C
ATOM	4838	O	ASP	A 324	-44.108	-9.885	6.255	1.00	33.36	O
ATOM	4840	N	TYR	A 325	-43.606	-9.423	8.392	1.00	32.10	N
ATOM	4841	CA	TYR	A 325	-44.515	-8.279	8.400	1.00	31.65	C
ATOM	4843	CB	TYR	A 325	-44.718	-7.726	9.815	1.00	31.68	C
ATOM	4846	CG	TYR	A 325	-43.618	-6.846	10.352	1.00	31.35	C
ATOM	4847	CD1	TYR	A 325	-42.507	-7.389	10.992	1.00	31.31	C
ATOM	4849	CE1	TYR	A 325	-41.497	-6.574	11.504	1.00	31.46	C
ATOM	4851	CZ	TYR	A 325	-41.613	-5.192	11.392	1.00	32.41	C
ATOM	4852	OH	TYR	A 325	-40.637	-4.336	11.893	1.00	32.78	O
ATOM	4854	CE2	TYR	A 325	-42.723	-4.644	10.769	1.00	32.17	C
ATOM	4856	CD2	TYR	A 325	-43.713	-5.470	10.261	1.00	31.47	C
ATOM	4858	C	TYR	A 325	-44.094	-7.179	7.424	1.00	31.17	C
ATOM	4859	O	TYR	A 325	-44.947	-6.458	6.920	1.00	31.09	O
ATOM	4861	N	MET	A 326	-42.796	-7.068	7.145	1.00	30.73	N
ATOM	4862	CA	MET	A 326	-42.277	-6.036	6.225	1.00	30.49	C
ATOM	4864	CB	MET	A 326	-40.832	-5.702	6.565	1.00	30.16	C
ATOM	4867	CG	MET	A 326	-40.725	-4.918	7.830	1.00	29.57	C
ATOM	4870	SD	MET	A 326	-39.057	-4.376	8.166	1.00	28.39	S
ATOM	4871	CE	MET	A 326	-39.407	-2.787	8.933	1.00	26.18	C
ATOM	4875	C	MET	A 326	-42.371	-6.418	4.748	1.00	30.59	C
ATOM	4876	O	MET	A 326	-42.786	-5.603	3.920	1.00	30.48	O
ATOM	4878	N	LYS	A 327	-41.936	-7.643	4.438	1.00	30.68	N
ATOM	4879	CA	LYS	A 327	-42.170	-8.299	3.143	1.00	30.56	C
ATOM	4881	CB	LYS	A 327	-42.267	-9.830	3.326	1.00	30.84	C
ATOM	4884	CG	LYS	A 327	-41.052	-10.607	2.844	1.00	32.35	C
ATOM	4887	CD	LYS	A 327	-40.990	-12.019	3.412	1.00	34.43	C
ATOM	4890	CE	LYS	A 327	-40.221	-12.941	2.462	1.00	35.77	C
ATOM	4893	NZ	LYS	A 327	-39.691	-14.162	3.150	1.00	37.20	N
ATOM	4897	C	LYS	A 327	-43.453	-7.811	2.515	1.00	29.89	C
ATOM	4898	O	LYS	A 327	-43.447	-7.220	1.436	1.00	29.85	O
ATOM	4900	N	LEU	A 328	-44.544	-8.046	3.230	1.00	29.10	N
ATOM	4901	CA	LEU	A 328	-45.871	-7.781	2.730	1.00	28.77	C
ATOM	4903	CB	LEU	A 328	-46.899	-8.336	3.709	1.00	28.81	C
ATOM	4906	CG	LEU	A 328	-48.349	-8.364	3.257	1.00	28.82	C
ATOM	4908	CD1	LEU	A 328	-48.519	-9.260	2.046	1.00	29.09	C
ATOM	4912	CD2	LEU	A 328	-49.202	-8.845	4.411	1.00	29.16	C
ATOM	4916	C	LEU	A 328	-46.057	-6.291	2.564	1.00	28.61	C
ATOM	4917	O	LEU	A 328	-46.582	-5.828	1.554	1.00	29.00	O
ATOM	4919	N	CYS	A 329	-45.612	-5.532	3.557	1.00	28.24	N
ATOM	4920	CA	CYS	A 329	-45.737	-4.086	3.504	1.00	27.95	C
ATOM	4922	CB	CYS	A 329	-45.311	-3.459	4.834	1.00	28.21	C
ATOM	4925	SG	CYS	A 329	-45.280	-1.630	4.817	1.00	32.39	S
ATOM	4927	C	CYS	A 329	-44.921	-3.541	2.327	1.00	26.34	C
ATOM	4928	O	CYS	A 329	-45.459	-2.853	1.475	1.00	26.23	O
ATOM	4930	N	PHE	A 330	-43.642	-3.882	2.269	1.00	24.82	N
ATOM	4931	CA	PHE	A 330	-42.790	-3.514	1.130	1.00	23.82	C
ATOM	4933	CB	PHE	A 330	-41.397	-4.137	1.283	1.00	23.64	C
ATOM	4936	CG	PHE	A 330	-40.492	-3.873	.117	1.00	22.63	C
ATOM	4937	CD1	PHE	A 330	-39.845	-2.658	-.008	1.00	21.68	C
ATOM	4939	CE1	PHE	A 330	-39.020	-2.393	-1.082	1.00	21.07	C
ATOM	4941	CZ	PHE	A 330	-38.829	-3.343	-2.046	1.00	22.00	C
ATOM	4943	CE2	PHE	A 330	-39.474	-4.568	-1.944	1.00	22.60	C
ATOM	4945	CD2	PHE	A 330	-40.309	-4.824	-.865	1.00	22.36	C
ATOM	4947	C	PHE	A 330	-43.348	-3.863	-.281	1.00	22.99	C
ATOM	4948	O	PHE	A 330	-43.350	-3.012	-1.182	1.00	22.98	O
ATOM	4950	N	LEU	A 331	-43.789	-5.101	-.492	1.00	21.65	N

ATOM	4951	CA	LEU A 331	-44.306	-5.484	-1.816	1.00	20.74	C
ATOM	4953	CB	LEU A 331	-44.573	-6.990	-1.912	1.00	20.57	C
ATOM	4956	CG	LEU A 331	-44.959	-7.575	-3.277	1.00	19.82	C
ATOM	4958	CD1	LEU A 331	-43.936	-7.246	-4.329	1.00	19.00	C
ATOM	4962	CD2	LEU A 331	-45.128	-9.092	-3.165	1.00	19.11	C
ATOM	4966	C	LEU A 331	-45.568	-4.710	-2.159	1.00	19.98	C
ATOM	4967	O	LEU A 331	-45.753	-4.324	-3.300	1.00	20.15	O
ATOM	4969	N	ALA A 332	-46.431	-4.495	-1.172	1.00	19.09	N
ATOM	4970	CA	ALA A 332	-47.619	-3.675	-1.353	1.00	18.46	C
ATOM	4972	CB	ALA A 332	-48.406	-3.587	-.045	1.00	18.20	C
ATOM	4976	C	ALA A 332	-47.248	-2.279	-1.856	1.00	17.91	C
ATOM	4977	O	ALA A 332	-47.890	-1.745	-2.744	1.00	17.98	O
ATOM	4979	N	LEU A 333	-46.197	-1.706	-1.295	1.00	17.62	N
ATOM	4980	CA	LEU A 333	-45.753	-.353	-1.637	1.00	17.60	C
ATOM	4982	CB	LEU A 333	-44.725	.132	-.598	1.00	17.59	C
ATOM	4985	CG	LEU A 333	-44.122	1.533	-.761	1.00	17.16	C
ATOM	4987	CD1	LEU A 333	-45.166	2.629	-.554	1.00	16.99	C
ATOM	4991	CD2	LEU A 333	-42.979	1.704	.200	1.00	15.89	C
ATOM	4995	C	LEU A 333	-45.100	-.320	-3.005	1.00	17.70	C
ATOM	4996	O	LEU A 333	-45.321	.603	-3.795	1.00	17.86	O
ATOM	4998	N	TYR A 334	-44.248	-1.319	-3.234	1.00	17.63	N
ATOM	4999	CA	TYR A 334	-43.542	-1.531	-4.489	1.00	17.44	C
ATOM	5001	CB	TYR A 334	-42.893	-2.908	-4.444	1.00	17.50	C
ATOM	5004	CG	TYR A 334	-41.897	-3.169	-5.523	1.00	18.03	C
ATOM	5005	CD1	TYR A 334	-40.698	-2.473	-5.569	1.00	18.83	C
ATOM	5007	CE1	TYR A 334	-39.764	-2.725	-6.542	1.00	18.13	C
ATOM	5009	CZ	TYR A 334	-40.016	-3.678	-7.478	1.00	18.78	C
ATOM	5010	OH	TYR A 334	-39.096	-3.929	-8.439	1.00	21.71	O
ATOM	5012	CE2	TYR A 334	-41.180	-4.391	-7.462	1.00	19.74	C
ATOM	5014	CD2	TYR A 334	-42.118	-4.139	-6.472	1.00	19.59	C
ATOM	5016	C	TYR A 334	-44.500	-1.497	-5.649	1.00	17.43	C
ATOM	5017	O	TYR A 334	-44.264	-.803	-6.636	1.00	17.64	O
ATOM	5019	N	ASN A 335	-45.589	-2.253	-5.504	1.00	17.30	N
ATOM	5020	CA	ASN A 335	-46.574	-2.448	-6.553	1.00	17.23	C
ATOM	5022	CB	ASN A 335	-47.544	-3.551	-6.166	1.00	17.32	C
ATOM	5025	CG	ASN A 335	-46.952	-4.920	-6.332	1.00	18.16	C
ATOM	5026	OD1	ASN A 335	-45.913	-5.090	-6.989	1.00	18.99	O
ATOM	5027	ND2	ASN A 335	-47.616	-5.921	-5.749	1.00	18.49	N
ATOM	5030	C	ASN A 335	-47.365	-1.218	-6.812	1.00	17.12	C
ATOM	5031	O	ASN A 335	-47.613	-.852	-7.965	1.00	17.36	O
ATOM	5033	N	THR A 336	-47.789	-.602	-5.722	1.00	17.20	N
ATOM	5034	CA	THR A 336	-48.601	.593	-5.779	1.00	17.36	C
ATOM	5036	CB	THR A 336	-48.888	1.104	-4.381	1.00	17.20	C
ATOM	5038	OG1	THR A 336	-49.611	.103	-3.657	1.00	16.26	O
ATOM	5040	CG2	THR A 336	-49.688	2.384	-4.452	1.00	17.18	C
ATOM	5044	C	THR A 336	-47.893	1.691	-6.550	1.00	17.84	C
ATOM	5045	O	THR A 336	-48.511	2.388	-7.360	1.00	18.07	O
ATOM	5047	N	ILE A 337	-46.595	1.831	-6.298	1.00	18.14	N
ATOM	5048	CA	ILE A 337	-45.817	2.904	-6.895	1.00	18.50	C
ATOM	5050	CB	ILE A 337	-44.584	3.209	-6.066	1.00	18.31	C
ATOM	5052	CG1	ILE A 337	-45.014	3.997	-4.837	1.00	18.58	C
ATOM	5055	CD1	ILE A 337	-44.043	3.866	-3.735	1.00	20.66	C
ATOM	5059	CG2	ILE A 337	-43.570	3.988	-6.867	1.00	16.98	C
ATOM	5063	C	ILE A 337	-45.447	2.548	-8.314	1.00	19.29	C
ATOM	5064	O	ILE A 337	-45.556	3.387	-9.214	1.00	19.17	O
ATOM	5066	N	ASN A 338	-45.033	1.299	-8.513	1.00	20.16	N
ATOM	5067	CA	ASN A 338	-44.861	.767	-9.864	1.00	20.93	C
ATOM	5069	CB	ASN A 338	-44.409	-.695	-9.830	1.00	21.07	C
ATOM	5072	CG	ASN A 338	-42.953	-.845	-9.439	1.00	21.61	C
ATOM	5073	OD1	ASN A 338	-42.232	.143	-9.308	1.00	22.40	O

ATOM	5074	ND2 ASN A 338	-42.509	-2.086	-9.260	1.00	21.92	N
ATOM	5077	C ASN A 338	-46.123	.914	-10.719	1.00	21.44	C
ATOM	5078	O ASN A 338	-46.022	1.202	-11.904	1.00	21.25	O
ATOM	5080	N GLU A 339	-47.303	.747	-10.128	1.00	22.21	N
ATOM	5081	CA GLU A 339	-48.532	1.000	-10.880	1.00	23.31	C
ATOM	5083	CB GLU A 339	-49.768	.464	-10.164	1.00	24.13	C
ATOM	5086	CG GLU A 339	-50.146	-.936	-10.660	1.00	28.69	C
ATOM	5089	CD GLU A 339	-50.939	-1.740	-9.638	1.00	35.06	C
ATOM	5090	OE1 GLU A 339	-51.899	-1.136	-9.067	1.00	39.51	O
ATOM	5091	OE2 GLU A 339	-50.596	-2.952	-9.414	1.00	36.64	O
ATOM	5092	C GLU A 339	-48.725	2.463	-11.269	1.00	22.71	C
ATOM	5093	O GLU A 339	-49.215	2.743	-12.372	1.00	22.48	O
ATOM	5095	N ILE A 340	-48.339	3.390	-10.390	1.00	22.19	N
ATOM	5096	CA ILE A 340	-48.406	4.813	-10.734	1.00	21.64	C
ATOM	5098	CB ILE A 340	-48.135	5.739	-9.538	1.00	21.33	C
ATOM	5100	CG1 ILE A 340	-49.229	5.596	-8.482	1.00	21.29	C
ATOM	5103	CD1 ILE A 340	-48.925	6.306	-7.150	1.00	20.00	C
ATOM	5107	CG2 ILE A 340	-48.091	7.177	-9.982	1.00	20.28	C
ATOM	5111	C ILE A 340	-47.414	5.128	-11.861	1.00	21.78	C
ATOM	5112	O ILE A 340	-47.786	5.826	-12.818	1.00	22.13	O
ATOM	5114	N ALA A 341	-46.179	4.609	-11.771	1.00	21.34	N
ATOM	5115	CA ALA A 341	-45.147	4.892	-12.787	1.00	21.03	C
ATOM	5117	CB ALA A 341	-43.837	4.268	-12.406	1.00	20.55	C
ATOM	5121	C ALA A 341	-45.592	4.426	-14.183	1.00	21.32	C
ATOM	5122	O ALA A 341	-45.228	5.036	-15.217	1.00	21.11	O
ATOM	5124	N TYR A 342	-46.393	3.355	-14.196	1.00	21.45	N
ATOM	5125	CA TYR A 342	-47.008	2.861	-15.414	1.00	21.41	C
ATOM	5127	CB TYR A 342	-47.627	1.468	-15.208	1.00	21.31	C
ATOM	5130	CG TYR A 342	-48.336	.957	-16.450	1.00	19.60	C
ATOM	5131	CD1 TYR A 342	-47.613	.463	-17.521	1.00	16.69	C
ATOM	5133	CE1 TYR A 342	-48.231	.032	-18.643	1.00	15.50	C
ATOM	5135	CZ TYR A 342	-49.595	.085	-18.728	1.00	15.88	C
ATOM	5136	OH TYR A 342	-50.196	-.352	-19.877	1.00	16.46	O
ATOM	5138	CE2 TYR A 342	-50.352	.558	-17.681	1.00	16.50	C
ATOM	5140	CD2 TYR A 342	-49.725	.997	-16.556	1.00	18.17	C
ATOM	5142	C TYR A 342	-48.064	3.830	-15.899	1.00	22.08	C
ATOM	5143	O TYR A 342	-48.094	4.169	-17.048	1.00	22.07	O
ATOM	5145	N ASP A 343	-48.942	4.283	-15.032	1.00	23.41	N
ATOM	5146	CA ASP A 343	-49.970	5.211	-15.481	1.00	24.66	C
ATOM	5148	CB ASP A 343	-50.851	5.674	-14.318	1.00	25.10	C
ATOM	5151	CG ASP A 343	-51.720	4.552	-13.732	1.00	26.06	C
ATOM	5152	OD1 ASP A 343	-52.107	3.613	-14.477	1.00	26.39	O
ATOM	5153	OD2 ASP A 343	-52.032	4.642	-12.516	1.00	27.29	O
ATOM	5154	C ASP A 343	-49.316	6.420	-16.142	1.00	25.22	C
ATOM	5155	O ASP A 343	-49.755	6.874	-17.192	1.00	25.47	O
ATOM	5157	N ASN A 344	-48.260	6.936	-15.528	1.00	25.87	N
ATOM	5158	CA ASN A 344	-47.508	8.042	-16.127	1.00	26.41	C
ATOM	5160	CB ASN A 344	-46.498	8.605	-15.134	1.00	26.56	C
ATOM	5163	CG ASN A 344	-47.152	9.407	-14.073	1.00	26.95	C
ATOM	5164	OD1 ASN A 344	-47.495	10.568	-14.296	1.00	29.13	O
ATOM	5165	ND2 ASN A 344	-47.367	8.799	-12.916	1.00	26.36	N
ATOM	5168	C ASN A 344	-46.785	7.674	-17.416	1.00	26.53	C
ATOM	5169	O ASN A 344	-46.658	8.508	-18.304	1.00	26.74	O
ATOM	5171	N LEU A 345	-46.280	6.448	-17.510	1.00	26.56	N
ATOM	5172	CA LEU A 345	-45.634	6.019	-18.739	1.00	26.37	C
ATOM	5174	CB LEU A 345	-44.890	4.693	-18.550	1.00	26.32	C
ATOM	5177	CG LEU A 345	-43.995	4.339	-19.750	1.00	25.96	C
ATOM	5179	CD1 LEU A 345	-42.706	5.133	-19.665	1.00	25.07	C
ATOM	5183	CD2 LEU A 345	-43.724	2.835	-19.874	1.00	24.75	C
ATOM	5187	C LEU A 345	-46.679	5.908	-19.856	1.00	26.38	C

ATOM	5188	O	LEU A 345	-46.435	6.335	-20.966	1.00	26.52	O
ATOM	5190	N	LYS A 346	-47.839	5.340	-19.563	1.00	26.59	N
ATOM	5191	CA	LYS A 346	-48.880	5.176	-20.572	1.00	26.89	C
ATOM	5193	CB	LYS A 346	-50.065	4.367	-20.011	1.00	26.80	C
ATOM	5196	CG	LYS A 346	-51.073	3.931	-21.062	1.00	26.26	C
ATOM	5199	CD	LYS A 346	-52.210	3.080	-20.517	1.00	26.15	C
ATOM	5202	CE	LYS A 346	-53.227	3.849	-19.689	1.00	26.40	C
ATOM	5205	NZ	LYS A 346	-53.136	3.506	-18.223	1.00	27.67	N
ATOM	5209	C	LYS A 346	-49.372	6.537	-21.071	1.00	27.52	C
ATOM	5210	O	LYS A 346	-49.562	6.742	-22.272	1.00	27.57	O
ATOM	5212	N	ASP A 347	-49.567	7.472	-20.148	1.00	28.10	N
ATOM	5213	CA	ASP A 347	-50.309	8.689	-20.465	1.00	28.58	C
ATOM	5215	CB	ASP A 347	-51.305	9.017	-19.329	1.00	28.85	C
ATOM	5218	CG	ASP A 347	-52.426	7.950	-19.197	1.00	30.01	C
ATOM	5219	OD1	ASP A 347	-52.827	7.349	-20.223	1.00	30.51	O
ATOM	5220	OD2	ASP A 347	-52.910	7.704	-18.069	1.00	32.61	O
ATOM	5221	C	ASP A 347	-49.407	9.871	-20.804	1.00	28.22	C
ATOM	5222	O	ASP A 347	-49.778	10.710	-21.611	1.00	28.39	O
ATOM	5224	N	LYS A 348	-48.228	9.930	-20.206	1.00	28.03	N
ATOM	5225	CA	LYS A 348	-47.301	11.021	-20.467	1.00	28.01	C
ATOM	5227	CB	LYS A 348	-46.785	11.648	-19.164	1.00	28.39	C
ATOM	5230	CG	LYS A 348	-47.834	12.305	-18.257	1.00	30.18	C
ATOM	5233	CD	LYS A 348	-47.143	12.898	-17.003	1.00	32.40	C
ATOM	5236	CE	LYS A 348	-48.124	13.121	-15.874	1.00	33.59	C
ATOM	5239	NZ	LYS A 348	-49.319	13.878	-16.341	1.00	35.48	N
ATOM	5243	C	LYS A 348	-46.107	10.551	-21.260	1.00	27.25	C
ATOM	5244	O	LYS A 348	-45.241	11.345	-21.566	1.00	27.60	O
ATOM	5246	N	GLY A 349	-46.036	9.271	-21.583	1.00	26.47	N
ATOM	5247	CA	GLY A 349	-44.863	8.742	-22.261	1.00	26.01	C
ATOM	5250	C	GLY A 349	-43.559	9.109	-21.587	1.00	25.67	C
ATOM	5251	O	GLY A 349	-42.613	9.472	-22.250	1.00	25.85	O
ATOM	5253	N	GLU A 350	-43.498	9.032	-20.269	1.00	25.55	N
ATOM	5254	CA	GLU A 350	-42.272	9.374	-19.555	1.00	25.76	C
ATOM	5256	CB	GLU A 350	-42.332	10.814	-19.021	1.00	26.34	C
ATOM	5259	CG	GLU A 350	-42.179	11.928	-20.106	1.00	29.04	C
ATOM	5262	CD	GLU A 350	-40.745	12.101	-20.598	1.00	32.08	C
ATOM	5263	OE1	GLU A 350	-39.853	12.170	-19.716	1.00	35.46	O
ATOM	5264	OE2	GLU A 350	-40.516	12.174	-21.839	1.00	31.64	O
ATOM	5265	C	GLU A 350	-42.084	8.414	-18.400	1.00	24.87	C
ATOM	5266	O	GLU A 350	-43.067	8.046	-17.760	1.00	25.16	O
ATOM	5268	N	ASN A 351	-40.833	8.009	-18.145	1.00	23.71	N
ATOM	5269	CA	ASN A 351	-40.494	7.189	-16.980	1.00	22.88	C
ATOM	5271	CB	ASN A 351	-39.351	6.229	-17.287	1.00	22.88	C
ATOM	5274	CG	ASN A 351	-39.141	5.215	-16.184	1.00	22.98	C
ATOM	5275	OD1	ASN A 351	-39.006	5.584	-15.035	1.00	22.68	O
ATOM	5276	ND1	ASN A 351	-39.139	3.924	-16.529	1.00	24.38	N
ATOM	5279	C	ASN A 351	-40.126	8.058	-15.786	1.00	22.40	C
ATOM	5280	O	ASN A 351	-39.097	8.741	-15.785	1.00	22.66	O
ATOM	5282	N	ILE A 352	-40.965	8.021	-14.760	1.00	21.55	N
ATOM	5283	CA	ILE A 352	-40.770	8.860	-13.591	1.00	20.75	C
ATOM	5285	CB	ILE A 352	-42.002	9.788	-13.360	1.00	20.95	C
ATOM	5287	CG1	ILE A 352	-43.263	8.970	-13.022	1.00	20.82	C
ATOM	5290	CD1	ILE A 352	-44.237	9.698	-12.122	1.00	20.21	C
ATOM	5294	CG2	ILE A 352	-42.266	10.669	-14.582	1.00	19.85	C
ATOM	5298	C	ILE A 352	-40.516	8.016	-12.339	1.00	20.37	C
ATOM	5299	O	ILE A 352	-40.574	8.513	-11.224	1.00	20.55	O
ATOM	5301	N	LEU A 353	-40.236	6.733	-12.519	1.00	19.98	N
ATOM	5302	CA	LEU A 353	-40.061	5.828	-11.387	1.00	19.42	C
ATOM	5304	CB	LEU A 353	-39.894	4.392	-11.870	1.00	19.27	C
ATOM	5307	CG	LEU A 353	-39.989	3.318	-10.799	1.00	19.01	C

ATOM	5309	CD1 LEU A 353	-41.164	3.574	-9.867	1.00	18.35	C
ATOM	5313	CD2 LEU A 353	-40.085	1.938	-11.468	1.00	18.74	C
ATOM	5317	C LEU A 353	-38.890	6.240	-10.523	1.00	19.11	C
ATOM	5318	O LEU A 353	-39.022	6.319	-9.314	1.00	18.96	O
ATOM	5320	N PRO A 354	-37.746	6.549	-11.142	1.00	19.06	N
ATOM	5321	CA PRO A 354	-36.631	7.069	-10.367	1.00	19.29	C
ATOM	5323	CB PRO A 354	-35.749	7.717	-11.438	1.00	19.25	C
ATOM	5326	CG PRO A 354	-35.997	6.942	-12.638	1.00	19.08	C
ATOM	5329	CD PRO A 354	-37.409	6.472	-12.572	1.00	19.00	C
ATOM	5332	C PRO A 354	-37.024	8.122	-9.328	1.00	19.33	C
ATOM	5333	O PRO A 354	-36.534	8.057	-8.193	1.00	19.73	O
ATOM	5334	N TYR A 355	-37.891	9.063	-9.721	1.00	18.83	N
ATOM	5335	CA TYR A 355	-38.206	10.221	-8.898	1.00	18.79	C
ATOM	5337	CB TYR A 355	-38.836	11.359	-9.707	1.00	19.10	C
ATOM	5340	CG TYR A 355	-38.142	11.708	-11.009	1.00	20.15	C
ATOM	5341	CD1 TYR A 355	-36.940	12.403	-11.027	1.00	20.78	C
ATOM	5343	CE1 TYR A 355	-36.321	12.719	-12.232	1.00	22.19	C
ATOM	5345	CZ TYR A 355	-36.918	12.347	-13.438	1.00	22.19	C
ATOM	5346	OH TYR A 355	-36.341	12.660	-14.657	1.00	23.29	O
ATOM	5348	CE2 TYR A 355	-38.113	11.677	-13.430	1.00	21.62	C
ATOM	5350	CD2 TYR A 355	-38.720	11.371	-12.227	1.00	21.05	C
ATOM	5352	C TYR A 355	-39.155	9.864	-7.777	1.00	18.43	C
ATOM	5353	O TYR A 355	-39.081	10.460	-6.709	1.00	18.61	O
ATOM	5355	N LEU A 356	-40.058	8.917	-8.023	1.00	17.92	N
ATOM	5356	CA LEU A 356	-40.993	8.467	-6.994	1.00	17.38	C
ATOM	5358	CB LEU A 356	-42.136	7.658	-7.597	1.00	17.17	C
ATOM	5361	CG LEU A 356	-42.956	8.366	-8.682	1.00	17.56	C
ATOM	5363	CD1 LEU A 356	-43.933	7.419	-9.371	1.00	17.47	C
ATOM	5367	CD2 LEU A 356	-43.698	9.548	-8.112	1.00	18.08	C
ATOM	5371	C LEU A 356	-40.241	7.623	-5.978	1.00	17.10	C
ATOM	5372	O LEU A 356	-40.332	7.859	-4.783	1.00	17.62	O
ATOM	5374	N THR A 357	-39.464	6.656	-6.442	1.00	16.64	N
ATOM	5375	CA THR A 357	-38.775	5.769	-5.513	1.00	16.18	C
ATOM	5377	CB THR A 357	-38.104	4.554	-6.213	1.00	16.04	C
ATOM	5379	OG1 THR A 357	-37.092	4.996	-7.123	1.00	15.93	O
ATOM	5381	CG2 THR A 357	-39.142	3.732	-6.962	1.00	14.99	C
ATOM	5385	C THR A 357	-37.764	6.544	-4.686	1.00	16.16	C
ATOM	5386	O THR A 357	-37.599	6.278	-3.506	1.00	16.07	O
ATOM	5388	N LYS A 358	-37.107	7.524	-5.291	1.00	16.20	N
ATOM	5389	CA LYS A 358	-36.186	8.365	-4.543	1.00	16.46	C
ATOM	5391	CB LYS A 358	-35.518	9.386	-5.453	1.00	16.80	C
ATOM	5394	CG LYS A 358	-34.612	10.388	-4.741	1.00	18.30	C
ATOM	5397	CD LYS A 358	-33.352	9.738	-4.168	1.00	20.21	C
ATOM	5400	CE LYS A 358	-32.335	10.811	-3.768	1.00	21.96	C
ATOM	5403	NZ LYS A 358	-31.163	10.275	-3.019	1.00	22.95	N
ATOM	5407	C LYS A 358	-36.936	9.083	-3.440	1.00	16.26	C
ATOM	5408	O LYS A 358	-36.448	9.161	-2.320	1.00	16.15	O
ATOM	5410	N ALA A 359	-38.126	9.593	-3.763	1.00	16.11	N
ATOM	5411	CA ALA A 359	-38.937	10.332	-2.798	1.00	16.10	C
ATOM	5413	CB ALA A 359	-40.221	10.768	-3.406	1.00	15.71	C
ATOM	5417	C ALA A 359	-39.215	9.476	-1.588	1.00	16.52	C
ATOM	5418	O ALA A 359	-39.247	9.970	-.442	1.00	16.82	O
ATOM	5420	N TRP A 360	-39.398	8.187	-1.843	1.00	16.79	N
ATOM	5421	CA TRP A 360	-39.704	7.247	-.780	1.00	17.25	C
ATOM	5423	CB TRP A 360	-40.390	6.006	-1.352	1.00	17.41	C
ATOM	5426	CG TRP A 360	-41.852	6.129	-1.318	1.00	17.06	C
ATOM	5427	CD1 TRP A 360	-42.664	6.451	-2.346	1.00	17.94	C
ATOM	5429	NE1 TRP A 360	-43.967	6.489	-1.926	1.00	18.03	N
ATOM	5431	CE2 TRP A 360	-44.002	6.193	-.592	1.00	17.72	C
ATOM	5432	CD2 TRP A 360	-42.684	5.965	-.179	1.00	16.60	C

ATOM	5433	CE3 TRP A 360	-42.441	5.654	1.155	1.00	16.65	C
ATOM	5435	CZ3 TRP A 360	-43.508	5.573	2.022	1.00	16.91	C
ATOM	5437	CH2 TRP A 360	-44.811	5.807	1.586	1.00	17.82	C
ATOM	5439	CZ2 TRP A 360	-45.080	6.114	.282	1.00	18.49	C
ATOM	5441	C TRP A 360	-38.490	6.865	.073	1.00	17.52	C
ATOM	5442	O TRP A 360	-38.603	6.767	1.297	1.00	17.38	O
ATOM	5444	N ALA A 361	-37.344	6.651	-.568	1.00	17.82	N
ATOM	5445	CA ALA A 361	-36.122	6.336	.157	1.00	18.11	C
ATOM	5447	CB ALA A 361	-34.982	6.050	-.805	1.00	17.92	C
ATOM	5451	C ALA A 361	-35.781	7.507	1.063	1.00	18.52	C
ATOM	5452	O ALA A 361	-35.434	7.327	2.229	1.00	18.48	O
ATOM	5454	N ASP A 362	-35.911	8.711	.521	1.00	19.16	N
ATOM	5455	CA ASP A 362	-35.627	9.925	1.276	1.00	19.90	C
ATOM	5457	CB ASP A 362	-35.797	11.167	.387	1.00	20.39	C
ATOM	5460	CG ASP A 362	-34.596	11.424	-.530	1.00	21.84	C
ATOM	5461	OD1 ASP A 362	-33.630	10.625	-.561	1.00	22.88	O
ATOM	5462	OD2 ASP A 362	-34.630	12.453	-1.231	1.00	24.83	O
ATOM	5463	C ASP A 362	-36.532	10.039	2.510	1.00	19.91	C
ATOM	5464	O ASP A 362	-36.074	10.427	3.591	1.00	19.76	O
ATOM	5466	N LEU A 363	-37.813	9.707	2.346	1.00	19.95	N
ATOM	5467	CA LEU A 363	-38.747	9.688	3.478	1.00	19.85	C
ATOM	5469	CB LEU A 363	-40.175	9.415	3.006	1.00	19.70	C
ATOM	5472	CG LEU A 363	-41.219	9.293	4.123	1.00	18.58	C
ATOM	5474	CD1 LEU A 363	-41.189	10.551	4.965	1.00	18.37	C
ATOM	5478	CD2 LEU A 363	-42.597	9.073	3.551	1.00	16.43	C
ATOM	5482	C LEU A 363	-38.368	8.613	4.488	1.00	20.24	C
ATOM	5483	O LEU A 363	-38.314	8.875	5.691	1.00	20.52	O
ATOM	5485	N CYS A 364	-38.129	7.397	3.997	1.00	20.29	N
ATOM	5486	CA CYS A 364	-37.741	6.308	4.874	1.00	20.39	C
ATOM	5488	CB CYS A 364	-37.595	4.985	4.111	1.00	20.37	C
ATOM	5491	SG CYS A 364	-39.208	4.201	3.666	1.00	20.69	S
ATOM	5493	C CYS A 364	-36.467	6.683	5.646	1.00	20.59	C
ATOM	5494	O CYS A 364	-36.386	6.445	6.863	1.00	20.82	O
ATOM	5496	N ASN A 365	-35.495	7.314	4.980	1.00	20.35	N
ATOM	5497	CA ASN A 365	-34.282	7.716	5.697	1.00	20.14	C
ATOM	5499	CB ASN A 365	-33.188	8.203	4.754	1.00	20.18	C
ATOM	5502	CG ASN A 365	-32.359	7.064	4.184	1.00	20.55	C
ATOM	5503	OD1 ASN A 365	-31.706	6.316	4.925	1.00	20.40	O
ATOM	5504	ND2 ASN A 365	-32.365	6.938	2.854	1.00	21.25	N
ATOM	5507	C ASN A 365	-34.590	8.746	6.779	1.00	19.92	C
ATOM	5508	O ASN A 365	-33.997	8.685	7.857	1.00	19.94	O
ATOM	5510	N ALA A 366	-35.531	9.658	6.507	1.00	19.58	N
ATOM	5511	CA ALA A 366	-36.036	10.578	7.537	1.00	19.39	C
ATOM	5513	CB ALA A 366	-37.083	11.507	6.971	1.00	18.79	C
ATOM	5517	C ALA A 366	-36.597	9.784	8.730	1.00	19.75	C
ATOM	5518	O ALA A 366	-36.215	10.049	9.891	1.00	19.53	O
ATOM	5520	N PHE A 367	-37.460	8.797	8.447	1.00	19.69	N
ATOM	5521	CA PHE A 367	-37.985	7.925	9.502	1.00	20.11	C
ATOM	5523	CB PHE A 367	-38.952	6.857	8.967	1.00	20.47	C
ATOM	5526	CG PHE A 367	-40.293	7.370	8.494	1.00	21.35	C
ATOM	5527	CD1 PHE A 367	-40.985	8.347	9.180	1.00	21.54	C
ATOM	5529	CE1 PHE A 367	-42.224	8.775	8.732	1.00	21.47	C
ATOM	5531	CZ PHE A 367	-42.798	8.213	7.612	1.00	21.87	C
ATOM	5533	CE2 PHE A 367	-42.135	7.230	6.924	1.00	22.60	C
ATOM	5535	CD2 PHE A 367	-40.894	6.798	7.373	1.00	22.90	C
ATOM	5537	C PHE A 367	-36.871	7.179	10.252	1.00	20.12	C
ATOM	5538	O PHE A 367	-36.940	7.009	11.476	1.00	19.88	O
ATOM	5540	N LEU A 368	-35.868	6.696	9.516	1.00	20.30	N
ATOM	5541	CA LEU A 368	-34.775	5.931	10.133	1.00	20.19	C
ATOM	5543	CB LEU A 368	-33.783	5.413	9.085	1.00	20.08	C

ATOM	5546	CG	LEU A 368	-32.743	4.363	9.514	1.00	19.30	C
ATOM	5548	CD1	LEU A 368	-33.384	3.227	10.260	1.00	18.66	C
ATOM	5552	CD2	LEU A 368	-31.968	3.807	8.319	1.00	18.16	C
ATOM	5556	C	LEU A 368	-34.063	6.812	11.128	1.00	20.41	C
ATOM	5557	O	LEU A 368	-33.842	6.407	12.257	1.00	20.01	O
ATOM	5559	N	GLN A 369	-33.751	8.036	10.711	1.00	20.89	N
ATOM	5560	CA	GLN A 369	-33.037	8.970	11.564	1.00	21.49	C
ATOM	5562	CB	GLN A 369	-32.782	10.280	10.832	1.00	21.58	C
ATOM	5565	CG	GLN A 369	-32.071	11.359	11.677	1.00	21.20	C
ATOM	5568	CD	GLN A 369	-30.639	11.006	11.976	1.00	20.04	C
ATOM	5569	OE1	GLN A 369	-30.282	10.650	13.108	1.00	19.10	O
ATOM	5570	NE2	GLN A 369	-29.803	11.095	10.956	1.00	18.75	N
ATOM	5573	C	GLN A 369	-33.763	9.265	12.870	1.00	22.21	C
ATOM	5574	O	GLN A 369	-33.122	9.343	13.909	1.00	22.41	O
ATOM	5576	N	GLU A 370	-35.080	9.455	12.822	1.00	23.05	N
ATOM	5577	CA	GLU A 370	-35.856	9.718	14.046	1.00	23.71	C
ATOM	5579	CB	GLU A 370	-37.329	10.041	13.726	1.00	24.04	C
ATOM	5582	CG	GLU A 370	-37.484	11.293	12.862	1.00	26.84	C
ATOM	5585	CD	GLU A 370	-38.897	11.910	12.834	1.00	30.36	C
ATOM	5586	OE1	GLU A 370	-39.886	11.161	12.586	1.00	31.90	O
ATOM	5587	OE2	GLU A 370	-38.992	13.164	13.014	1.00	31.62	O
ATOM	5588	C	GLU A 370	-35.755	8.523	14.994	1.00	23.60	C
ATOM	5589	O	GLU A 370	-35.534	8.689	16.199	1.00	23.33	O
ATOM	5591	N	ALA A 371	-35.904	7.322	14.435	1.00	23.72	N
ATOM	5592	CA	ALA A 371	-35.771	6.091	15.201	1.00	23.91	C
ATOM	5594	CB	ALA A 371	-35.991	4.868	14.312	1.00	23.76	C
ATOM	5598	C	ALA A 371	-34.392	6.051	15.840	1.00	24.20	C
ATOM	5599	O	ALA A 371	-34.277	5.804	17.035	1.00	24.54	O
ATOM	5601	N	LYS A 372	-33.355	6.333	15.049	1.00	24.38	N
ATOM	5602	CA	LYS A 372	-31.979	6.273	15.530	1.00	24.39	C
ATOM	5604	CB	LYS A 372	-30.970	6.466	14.393	1.00	24.43	C
ATOM	5607	CG	LYS A 372	-30.623	5.164	13.645	1.00	25.20	C
ATOM	5610	CD	LYS A 372	-29.188	5.149	13.069	1.00	26.09	C
ATOM	5613	CE	LYS A 372	-29.114	5.437	11.557	1.00	27.01	C
ATOM	5616	NZ	LYS A 372	-29.002	4.200	10.709	1.00	26.82	N
ATOM	5620	C	LYS A 372	-31.717	7.268	16.645	1.00	24.59	C
ATOM	5621	O	LYS A 372	-31.096	6.908	17.627	1.00	25.13	O
ATOM	5623	N	TRP A 373	-32.181	8.507	16.520	1.00	24.78	N
ATOM	5624	CA	TRP A 373	-32.006	9.473	17.610	1.00	24.90	C
ATOM	5626	CB	TRP A 373	-32.565	10.863	17.266	1.00	24.75	C
ATOM	5629	CG	TRP A 373	-31.701	11.677	16.338	1.00	23.59	C
ATOM	5630	CD1	TRP A 373	-30.344	11.663	16.259	1.00	22.38	C
ATOM	5632	NE1	TRP A 373	-29.915	12.539	15.298	1.00	21.59	N
ATOM	5634	CE2	TRP A 373	-30.999	13.162	14.744	1.00	21.69	C
ATOM	5635	CD2	TRP A 373	-32.147	12.644	15.376	1.00	22.41	C
ATOM	5636	CE3	TRP A 373	-33.409	13.110	14.976	1.00	21.98	C
ATOM	5638	CZ3	TRP A 373	-33.480	14.072	13.979	1.00	22.08	C
ATOM	5640	CH2	TRP A 373	-32.313	14.567	13.366	1.00	22.31	C
ATOM	5642	CZ2	TRP A 373	-31.067	14.124	13.736	1.00	21.87	C
ATOM	5644	C	TRP A 373	-32.678	8.970	18.881	1.00	25.53	C
ATOM	5645	O	TRP A 373	-32.101	9.072	19.972	1.00	25.99	O
ATOM	5647	N	LEU A 374	-33.881	8.420	18.739	1.00	25.84	N
ATOM	5648	CA	LEU A 374	-34.667	7.980	19.893	1.00	26.29	C
ATOM	5650	CB	LEU A 374	-36.070	7.562	19.443	1.00	26.32	C
ATOM	5653	CG	LEU A 374	-37.227	7.646	20.444	1.00	26.39	C
ATOM	5655	CD1	LEU A 374	-38.456	8.286	19.761	1.00	26.95	C
ATOM	5659	CD2	LEU A 374	-37.583	6.281	21.050	1.00	26.24	C
ATOM	5663	C	LEU A 374	-33.982	6.824	20.623	1.00	26.77	C
ATOM	5664	O	LEU A 374	-33.995	6.758	21.860	1.00	26.95	O
ATOM	5666	N	TYR A 375	-33.383	5.919	19.854	1.00	27.23	N

ATOM	5667	CA	TYR A 375	-32.706	4.761	20.425	1.00	27.74	C
ATOM	5669	CB	TYR A 375	-32.195	3.814	19.328	1.00	27.78	C
ATOM	5672	CG	TYR A 375	-31.526	2.556	19.848	1.00	28.74	C
ATOM	5673	CD1	TYR A 375	-32.264	1.394	20.090	1.00	29.61	C
ATOM	5675	CE1	TYR A 375	-31.652	.232	20.573	1.00	29.89	C
ATOM	5677	CZ	TYR A 375	-30.287	.226	20.817	1.00	30.15	C
ATOM	5678	OH	TYR A 375	-29.676	-.917	21.286	1.00	30.57	O
ATOM	5680	CE2	TYR A 375	-29.530	1.367	20.579	1.00	30.05	C
ATOM	5682	CD2	TYR A 375	-30.152	2.522	20.096	1.00	29.74	C
ATOM	5684	C	TYR A 375	-31.553	5.248	21.275	1.00	27.97	C
ATOM	5685	O	TYR A 375	-31.404	4.837	22.422	1.00	28.00	O
ATOM	5687	N	ASN A 376	-30.763	6.158	20.719	1.00	28.32	N
ATOM	5688	CA	ASN A 376	-29.531	6.590	21.368	1.00	28.71	C
ATOM	5690	CB	ASN A 376	-28.569	7.161	20.329	1.00	28.60	C
ATOM	5693	CG	ASN A 376	-28.215	6.159	19.255	1.00	28.32	C
ATOM	5694	OD1	ASN A 376	-27.961	4.977	19.527	1.00	25.72	O
ATOM	5695	ND2	ASN A 376	-28.195	6.631	18.015	1.00	29.29	N
ATOM	5698	C	ASN A 376	-29.728	7.617	22.484	1.00	29.05	C
ATOM	5699	O	ASN A 376	-28.752	8.021	23.136	1.00	29.03	O
ATOM	5701	N	LYS A 377	-30.977	8.021	22.716	1.00	29.16	N
ATOM	5702	CA	LYS A 377	-31.254	9.178	23.549	1.00	29.27	C
ATOM	5704	CB	LYS A 377	-31.051	8.879	25.050	1.00	29.52	C
ATOM	5707	CG	LYS A 377	-32.202	8.112	25.723	1.00	30.36	C
ATOM	5710	CD	LYS A 377	-32.202	6.638	25.340	1.00	31.54	C
ATOM	5713	CE	LYS A 377	-33.322	5.862	26.031	1.00	32.30	C
ATOM	5716	NZ	LYS A 377	-33.559	4.509	25.411	1.00	32.13	N
ATOM	5720	C	LYS A 377	-30.337	10.299	23.080	1.00	28.92	C
ATOM	5721	O	LYS A 377	-29.590	10.867	23.875	1.00	28.94	O
ATOM	5723	N	SER A 378	-30.377	10.575	21.776	1.00	28.52	N
ATOM	5724	CA	SER A 378	-29.652	11.700	21.201	1.00	28.30	C
ATOM	5726	CB	SER A 378	-29.623	11.620	19.678	1.00	28.33	C
ATOM	5729	OG	SER A 378	-28.919	10.482	19.237	1.00	29.36	O
ATOM	5731	C	SER A 378	-30.355	12.977	21.594	1.00	27.97	C
ATOM	5732	O	SER A 378	-31.483	12.947	22.108	1.00	27.91	O
ATOM	5734	N	THR A 379	-29.684	14.098	21.344	1.00	27.60	N
ATOM	5735	CA	THR A 379	-30.257	15.422	21.581	1.00	27.24	C
ATOM	5737	CB	THR A 379	-29.929	15.950	23.002	1.00	27.23	C
ATOM	5739	OG1	THR A 379	-28.512	16.102	23.158	1.00	26.85	O
ATOM	5741	CG2	THR A 379	-30.467	15.007	24.067	1.00	27.48	C
ATOM	5745	C	THR A 379	-29.738	16.414	20.548	1.00	26.83	C
ATOM	5746	O	THR A 379	-28.834	17.190	20.844	1.00	26.69	O
ATOM	5748	N	PRO A 380	-30.305	16.390	19.331	1.00	26.56	N
ATOM	5749	CA	PRO A 380	-29.884	17.315	18.278	1.00	26.52	C
ATOM	5751	CB	PRO A 380	-30.380	16.648	16.992	1.00	26.46	C
ATOM	5754	CG	PRO A 380	-31.281	15.533	17.408	1.00	26.41	C
ATOM	5757	CD	PRO A 380	-31.395	15.508	18.887	1.00	26.53	C
ATOM	5760	C	PRO A 380	-30.469	18.728	18.386	1.00	26.55	C
ATOM	5761	O	PRO A 380	-31.472	18.966	19.063	1.00	26.69	O
ATOM	5762	N	THR A 381	-29.840	19.661	17.692	1.00	26.51	N
ATOM	5763	CA	THR A 381	-30.299	21.034	17.699	1.00	26.51	C
ATOM	5765	CB	THR A 381	-29.261	21.955	17.033	1.00	26.99	C
ATOM	5767	OG1	THR A 381	-28.919	21.431	15.732	1.00	27.65	O
ATOM	5769	CG2	THR A 381	-28.003	22.082	17.930	1.00	26.41	C
ATOM	5773	C	THR A 381	-31.635	21.167	16.970	1.00	26.06	C
ATOM	5774	O	THR A 381	-31.972	20.347	16.112	1.00	26.03	O
ATOM	5776	N	PHE A 382	-32.386	22.218	17.294	1.00	25.52	N
ATOM	5777	CA	PHE A 382	-33.681	22.436	16.654	1.00	24.83	C
ATOM	5779	CB	PHE A 382	-34.284	23.793	17.011	1.00	24.53	C
ATOM	5782	CG	PHE A 382	-35.495	24.113	16.211	1.00	23.77	C
ATOM	5783	CD1	PHE A 382	-36.745	23.715	16.637	1.00	24.50	C

ATOM	5785	CE1 PHE A 382	-37.873	23.977	15.878	1.00	24.49	C
ATOM	5787	CZ PHE A 382	-37.744	24.630	14.673	1.00	24.41	C
ATOM	5789	CE2 PHE A 382	-36.490	25.014	14.235	1.00	23.88	C
ATOM	5791	CD2 PHE A 382	-35.381	24.752	14.999	1.00	23.29	C
ATOM	5793	C PHE A 382	-33.542	22.342	15.150	1.00	24.49	C
ATOM	5794	O PHE A 382	-34.361	21.731	14.485	1.00	24.32	O
ATOM	5796	N ASP A 383	-32.498	22.969	14.627	1.00	24.42	N
ATOM	5797	CA ASP A 383	-32.257	23.005	13.190	1.00	24.19	C
ATOM	5799	CB ASP A 383	-31.101	23.966	12.867	1.00	24.21	C
ATOM	5802	CG ASP A 383	-31.473	25.423	13.050	1.00	23.72	C
ATOM	5803	OD1 ASP A 383	-32.634	25.794	12.837	1.00	25.60	O
ATOM	5804	OD2 ASP A 383	-30.594	26.220	13.383	1.00	24.03	O
ATOM	5805	C ASP A 383	-31.982	21.612	12.599	1.00	24.06	C
ATOM	5806	O ASP A 383	-32.393	21.342	11.473	1.00	24.00	O
ATOM	5808	N ASP A 384	-31.298	20.739	13.336	1.00	23.78	N
ATOM	5809	CA ASP A 384	-31.078	19.379	12.848	1.00	24.11	C
ATOM	5811	CB ASP A 384	-29.981	18.651	13.632	1.00	24.54	C
ATOM	5814	CG ASP A 384	-28.573	19.009	13.159	1.00	26.32	C
ATOM	5815	OD1 ASP A 384	-28.441	19.852	12.232	1.00	28.55	O
ATOM	5816	OD2 ASP A 384	-27.599	18.451	13.728	1.00	27.18	O
ATOM	5817	C ASP A 384	-32.351	18.550	12.905	1.00	23.83	C
ATOM	5818	O ASP A 384	-32.670	17.825	11.955	1.00	24.24	O
ATOM	5820	N TYR A 385	-33.070	18.637	14.017	1.00	23.31	N
ATOM	5821	CA TYR A 385	-34.294	17.861	14.179	1.00	22.89	C
ATOM	5823	CB TYR A 385	-34.833	17.979	15.608	1.00	22.86	C
ATOM	5826	CG TYR A 385	-36.144	17.245	15.807	1.00	22.65	C
ATOM	5827	CD1 TYR A 385	-36.163	15.868	16.030	1.00	22.36	C
ATOM	5829	CE1 TYR A 385	-37.351	15.189	16.206	1.00	22.77	C
ATOM	5831	CZ TYR A 385	-38.545	15.887	16.146	1.00	23.54	C
ATOM	5832	OH TYR A 385	-39.733	15.212	16.309	1.00	24.38	O
ATOM	5834	CE2 TYR A 385	-38.554	17.258	15.913	1.00	22.76	C
ATOM	5836	CD2 TYR A 385	-37.360	17.923	15.744	1.00	22.12	C
ATOM	5838	C TYR A 385	-35.375	18.299	13.190	1.00	22.55	C
ATOM	5839	O TYR A 385	-36.050	17.467	12.584	1.00	22.85	O
ATOM	5841	N PHE A 386	-35.537	19.609	13.042	1.00	22.08	N
ATOM	5842	CA PHE A 386	-36.626	20.155	12.248	1.00	21.62	C
ATOM	5844	CB PHE A 386	-36.857	21.628	12.568	1.00	21.73	C
ATOM	5847	CG PHE A 386	-38.033	22.209	11.851	1.00	21.54	C
ATOM	5848	CD1 PHE A 386	-39.311	21.933	12.273	1.00	21.29	C
ATOM	5850	CE1 PHE A 386	-40.380	22.444	11.623	1.00	22.11	C
ATOM	5852	CZ PHE A 386	-40.193	23.240	10.517	1.00	23.28	C
ATOM	5854	CE2 PHE A 386	-38.920	23.502	10.069	1.00	23.10	C
ATOM	5856	CD2 PHE A 386	-37.853	22.989	10.737	1.00	22.15	C
ATOM	5858	C PHE A 386	-36.369	19.987	10.769	1.00	21.09	C
ATOM	5859	O PHE A 386	-37.278	19.700	10.006	1.00	20.91	O
ATOM	5861	N GLY A 387	-35.127	20.177	10.362	1.00	20.68	N
ATOM	5862	CA GLY A 387	-34.747	19.917	8.985	1.00	20.43	C
ATOM	5865	C GLY A 387	-35.151	18.521	8.544	1.00	19.95	C
ATOM	5866	O GLY A 387	-35.553	18.327	7.398	1.00	20.24	O
ATOM	5868	N ASN A 388	-35.032	17.554	9.451	1.00	19.12	N
ATOM	5869	CA ASN A 388	-35.451	16.182	9.202	1.00	18.73	C
ATOM	5871	CB ASN A 388	-34.744	15.272	10.205	1.00	18.90	C
ATOM	5874	CG ASN A 388	-34.863	13.795	9.871	1.00	18.35	C
ATOM	5875	OD1 ASN A 388	-34.163	13.284	8.993	1.00	17.24	O
ATOM	5876	ND2 ASN A 388	-35.714	13.092	10.614	1.00	17.01	N
ATOM	5879	C ASN A 388	-36.968	16.033	9.350	1.00	18.65	C
ATOM	5880	O ASN A 388	-37.624	15.370	8.551	1.00	18.82	O
ATOM	5882	N ALA A 389	-37.527	16.670	10.371	1.00	18.38	N
ATOM	5883	CA ALA A 389	-38.926	16.477	10.728	1.00	18.02	C
ATOM	5885	CB ALA A 389	-39.203	17.107	12.084	1.00	17.83	C

ATOM	5889	C	ALA A 389	-39.949	16.966	9.700	1.00	17.90	C
ATOM	5890	O	ALA A 389	-41.095	16.544	9.760	1.00	17.74	O
ATOM	5892	N	TRP A 390	-39.585	17.860	8.783	1.00	18.11	N
ATOM	5893	CA	TRP A 390	-40.543	18.254	7.731	1.00	18.46	C
ATOM	5895	CB	TRP A 390	-40.459	19.737	7.329	1.00	18.59	C
ATOM	5898	CG	TRP A 390	-39.143	20.224	6.821	1.00	19.26	C
ATOM	5899	CD1	TRP A 390	-38.224	20.947	7.519	1.00	20.64	C
ATOM	5901	NE1	TRP A 390	-37.138	21.224	6.729	1.00	20.62	N
ATOM	5903	CE2	TRP A 390	-37.352	20.694	5.485	1.00	19.57	C
ATOM	5904	CD2	TRP A 390	-38.607	20.066	5.504	1.00	19.15	C
ATOM	5905	CE3	TRP A 390	-39.060	19.448	4.344	1.00	20.13	C
ATOM	5907	CZ3	TRP A 390	-38.254	19.477	3.224	1.00	20.27	C
ATOM	5909	CH2	TRP A 390	-37.017	20.110	3.239	1.00	19.38	C
ATOM	5911	CZ2	TRP A 390	-36.549	20.724	4.356	1.00	19.49	C
ATOM	5913	C	TRP A 390	-40.407	17.343	6.527	1.00	18.73	C
ATOM	5914	O	TRP A 390	-41.369	17.143	5.788	1.00	18.32	O
ATOM	5916	N	LYS A 391	-39.200	16.806	6.336	1.00	19.30	N
ATOM	5917	CA	LYS A 391	-38.969	15.692	5.404	1.00	19.70	C
ATOM	5919	CB	LYS A 391	-37.459	15.392	5.244	1.00	19.92	C
ATOM	5922	CG	LYS A 391	-36.812	15.953	3.935	1.00	22.08	C
ATOM	5925	CD	LYS A 391	-35.257	16.217	4.040	1.00	23.92	C
ATOM	5928	CE	LYS A 391	-34.918	17.698	4.426	1.00	24.43	C
ATOM	5931	NZ	LYS A 391	-33.699	17.901	5.307	1.00	23.23	N
ATOM	5935	C	LYS A 391	-39.731	14.441	5.861	1.00	19.46	C
ATOM	5936	O	LYS A 391	-40.284	13.709	5.032	1.00	20.07	O
ATOM	5938	N	SER A 392	-39.793	14.202	7.170	1.00	18.95	N
ATOM	5939	CA	SER A 392	-40.455	12.998	7.666	1.00	18.63	C
ATOM	5941	CB	SER A 392	-39.850	12.508	9.000	1.00	18.75	C
ATOM	5944	OG	SER A 392	-40.314	13.230	10.126	1.00	18.94	O
ATOM	5946	C	SER A 392	-41.964	13.152	7.771	1.00	18.20	C
ATOM	5947	O	SER A 392	-42.654	12.182	8.031	1.00	18.12	O
ATOM	5949	N	SER A 393	-42.477	14.360	7.567	1.00	18.07	N
ATOM	5950	CA	SER A 393	-43.929	14.591	7.556	1.00	18.02	C
ATOM	5952	CB	SER A 393	-44.229	16.078	7.504	1.00	18.08	C
ATOM	5955	OG	SER A 393	-43.995	16.558	6.192	1.00	18.01	O
ATOM	5957	C	SER A 393	-44.594	13.971	6.340	1.00	17.87	C
ATOM	5958	O	SER A 393	-45.778	13.661	6.372	1.00	18.03	O
ATOM	5960	N	SER A 394	-43.823	13.841	5.264	1.00	17.69	N
ATOM	5961	CA	SER A 394	-44.284	13.306	3.989	1.00	17.59	C
ATOM	5963	CB	SER A 394	-45.329	12.180	4.149	1.00	17.61	C
ATOM	5966	OG	SER A 394	-46.648	12.681	4.294	1.00	17.18	O
ATOM	5968	C	SER A 394	-44.828	14.415	3.115	1.00	17.42	C
ATOM	5969	O	SER A 394	-45.345	14.146	2.024	1.00	17.31	O
ATOM	5971	N	GLY A 395	-44.711	15.654	3.592	1.00	17.23	N
ATOM	5972	CA	GLY A 395	-45.088	16.827	2.807	1.00	17.26	C
ATOM	5975	C	GLY A 395	-44.478	16.749	1.415	1.00	17.29	C
ATOM	5976	O	GLY A 395	-45.203	16.653	.413	1.00	18.06	O
ATOM	5978	N	PRO A 396	-43.145	16.754	1.331	1.00	16.72	N
ATOM	5979	CA	PRO A 396	-42.582	16.610	.006	1.00	16.40	C
ATOM	5981	CB	PRO A 396	-41.080	16.549	.259	1.00	16.62	C
ATOM	5984	CG	PRO A 396	-40.903	16.749	1.775	1.00	17.11	C
ATOM	5987	CD	PRO A 396	-42.157	17.252	2.296	1.00	16.87	C
ATOM	5990	C	PRO A 396	-43.053	15.374	-.748	1.00	15.89	C
ATOM	5991	O	PRO A 396	-43.501	15.498	-1.894	1.00	15.90	O
ATOM	5992	N	LEU A 397	-42.973	14.197	-.135	1.00	15.27	N
ATOM	5993	CA	LEU A 397	-43.287	12.976	-.886	1.00	14.57	C
ATOM	5995	CB	LEU A 397	-43.332	11.733	-.008	1.00	14.36	C
ATOM	5998	CG	LEU A 397	-43.541	10.431	-.781	1.00	14.02	C
ATOM	6000	CD1	LEU A 397	-42.690	9.348	-.206	1.00	15.04	C
ATOM	6004	CD2	LEU A 397	-44.985	9.976	-.805	1.00	13.75	C

ATOM	6008	C	LEU A 397	-44.618	13.182	-1.542	1.00	14.39	C
ATOM	6009	O	LEU A 397	-44.736	12.999	-2.745	1.00	14.57	O
ATOM	6011	N	GLN A 398	-45.604	13.607	-.751	1.00	14.08	N
ATOM	6012	CA	GLN A 398	-46.962	13.828	-1.245	1.00	13.91	C
ATOM	6014	CB	GLN A 398	-47.860	14.363	-.136	1.00	13.99	C
ATOM	6017	CG	GLN A 398	-48.274	13.321	.909	1.00	14.01	C
ATOM	6020	CD	GLN A 398	-49.189	13.902	1.983	1.00	13.41	C
ATOM	6021	OE1	GLN A 398	-49.941	14.847	1.740	1.00	14.65	O
ATOM	6022	NE2	GLN A 398	-49.122	13.344	3.168	1.00	12.17	N
ATOM	6025	C	GLN A 398	-47.015	14.800	-2.403	1.00	13.84	C
ATOM	6026	O	GLN A 398	-47.677	14.547	-3.396	1.00	13.46	O
ATOM	6028	N	LEU A 399	-46.319	15.922	-2.272	1.00	14.11	N
ATOM	6029	CA	LEU A 399	-46.359	16.942	-3.315	1.00	14.37	C
ATOM	6031	CB	LEU A 399	-45.900	18.282	-2.756	1.00	14.26	C
ATOM	6034	CG	LEU A 399	-46.882	18.830	-1.704	1.00	14.44	C
ATOM	6036	CD1	LEU A 399	-46.250	19.960	-.919	1.00	16.43	C
ATOM	6040	CD2	LEU A 399	-48.191	19.304	-2.317	1.00	12.20	C
ATOM	6044	C	LEU A 399	-45.582	16.525	-4.578	1.00	14.68	C
ATOM	6045	O	LEU A 399	-46.043	16.784	-5.687	1.00	14.70	O
ATOM	6047	N	ILE A 400	-44.443	15.844	-4.419	1.00	14.92	N
ATOM	6048	CA	ILE A 400	-43.748	15.221	-5.564	1.00	15.08	C
ATOM	6050	CB	ILE A 400	-42.549	14.355	-5.129	1.00	15.35	C
ATOM	6052	CG1	ILE A 400	-41.406	15.254	-4.611	1.00	16.83	C
ATOM	6055	CD1	ILE A 400	-40.234	14.490	-3.958	1.00	17.18	C
ATOM	6059	CG2	ILE A 400	-42.060	13.486	-6.292	1.00	14.03	C
ATOM	6063	C	ILE A 400	-44.682	14.329	-6.358	1.00	15.08	C
ATOM	6064	O	ILE A 400	-44.672	14.362	-7.574	1.00	15.18	O
ATOM	6066	N	PHE A 401	-45.474	13.524	-5.650	1.00	15.16	N
ATOM	6067	CA	PHE A 401	-46.501	12.674	-6.258	1.00	14.90	C
ATOM	6069	CB	PHE A 401	-47.052	11.663	-5.240	1.00	14.51	C
ATOM	6072	CG	PHE A 401	-46.294	10.378	-5.201	1.00	12.92	C
ATOM	6073	CD1	PHE A 401	-46.727	9.282	-5.931	1.00	11.98	C
ATOM	6075	CE1	PHE A 401	-46.032	8.093	-5.917	1.00	11.20	C
ATOM	6077	CZ	PHE A 401	-44.881	7.985	-5.165	1.00	12.03	C
ATOM	6079	CE2	PHE A 401	-44.431	9.078	-4.428	1.00	11.98	C
ATOM	6081	CD2	PHE A 401	-45.143	10.263	-4.453	1.00	12.09	C
ATOM	6083	C	PHE A 401	-47.641	13.509	-6.806	1.00	15.39	C
ATOM	6084	O	PHE A 401	-48.183	13.215	-7.858	1.00	15.51	O
ATOM	6086	N	ALA A 402	-48.022	14.545	-6.080	1.00	16.15	N
ATOM	6087	CA	ALA A 402	-49.110	15.393	-6.525	1.00	16.90	C
ATOM	6089	CB	ALA A 402	-49.391	16.473	-5.505	1.00	16.99	C
ATOM	6093	C	ALA A 402	-48.722	16.003	-7.856	1.00	17.50	C
ATOM	6094	O	ALA A 402	-49.549	16.100	-8.770	1.00	17.48	O
ATOM	6096	N	TYR A 403	-47.444	16.367	-7.963	1.00	18.28	N
ATOM	6097	CA	TYR A 403	-46.916	17.046	-9.142	1.00	18.94	C
ATOM	6099	CB	TYR A 403	-45.412	17.252	-9.043	1.00	19.03	C
ATOM	6102	CG	TYR A 403	-44.823	17.801	-10.314	1.00	19.79	C
ATOM	6103	CD1	TYR A 403	-44.973	19.138	-10.652	1.00	20.80	C
ATOM	6105	CE1	TYR A 403	-44.436	19.643	-11.827	1.00	20.73	C
ATOM	6107	CZ	TYR A 403	-43.759	18.804	-12.674	1.00	21.16	C
ATOM	6108	OH	TYR A 403	-43.231	19.283	-13.836	1.00	22.98	O
ATOM	6110	CE2	TYR A 403	-43.608	17.476	-12.370	1.00	21.23	C
ATOM	6112	CD2	TYR A 403	-44.137	16.980	-11.195	1.00	20.93	C
ATOM	6114	C	TYR A 403	-47.198	16.293	-10.413	1.00	19.39	C
ATOM	6115	O	TYR A 403	-47.567	16.904	-11.422	1.00	19.57	O
ATOM	6117	N	PHE A 404	-47.023	14.976	-10.376	1.00	19.74	N
ATOM	6118	CA	PHE A 404	-47.239	14.177	-11.573	1.00	20.35	C
ATOM	6120	CB	PHE A 404	-46.533	12.841	-11.466	1.00	19.96	C
ATOM	6123	CG	PHE A 404	-45.048	12.971	-11.387	1.00	18.81	C
ATOM	6124	CD1	PHE A 404	-44.292	13.115	-12.528	1.00	17.42	C

ATOM	6126	CE1 PHE A 404	-42.935	13.239	-12.457	1.00	17.24	C
ATOM	6128	CZ PHE A 404	-42.312	13.237	-11.237	1.00	17.40	C
ATOM	6130	CE2 PHE A 404	-43.056	13.106	-10.093	1.00	17.72	C
ATOM	6132	CD2 PHE A 404	-44.413	12.975	-10.170	1.00	17.93	C
ATOM	6134	C PHE A 404	-48.713	13.994	-11.894	1.00	21.59	C
ATOM	6135	O PHE A 404	-49.072	13.717	-13.046	1.00	21.71	O
ATOM	6137	N ALA A 405	-49.572	14.169	-10.896	1.00	22.87	N
ATOM	6138	CA ALA A 405	-50.990	13.939	-11.099	1.00	24.06	C
ATOM	6140	CB ALA A 405	-51.604	13.408	-9.833	1.00	24.15	C
ATOM	6144	C ALA A 405	-51.724	15.192	-11.556	1.00	25.21	C
ATOM	6145	O ALA A 405	-52.876	15.114	-11.939	1.00	25.51	O
ATOM	6147	N VAL A 406	-51.056	16.336	-11.530	1.00	26.52	N
ATOM	6148	CA VAL A 406	-51.713	17.617	-11.759	1.00	27.68	C
ATOM	6150	CB VAL A 406	-51.694	18.406	-10.434	1.00	27.60	C
ATOM	6152	CG1 VAL A 406	-51.654	19.913	-10.663	1.00	28.07	C
ATOM	6156	CG2 VAL A 406	-52.883	18.009	-9.596	1.00	27.41	C
ATOM	6160	C VAL A 406	-51.097	18.429	-12.925	1.00	29.10	C
ATOM	6161	O VAL A 406	-51.810	19.104	-13.678	1.00	28.53	O
ATOM	6163	N VAL A 407	-49.772	18.356	-13.059	1.00	30.90	N
ATOM	6164	CA VAL A 407	-49.052	19.003	-14.151	1.00	32.15	C
ATOM	6166	CB VAL A 407	-47.574	19.228	-13.777	1.00	32.25	C
ATOM	6168	CG1 VAL A 407	-46.770	19.780	-14.964	1.00	32.22	C
ATOM	6172	CG2 VAL A 407	-47.490	20.154	-12.583	1.00	32.25	C
ATOM	6176	C VAL A 407	-49.134	18.160	-15.421	1.00	33.36	C
ATOM	6177	O VAL A 407	-48.688	17.009	-15.454	1.00	33.40	O
ATOM	6179	N GLN A 408	-49.693	18.767	-16.463	1.00	34.93	N
ATOM	6180	CA GLN A 408	-49.930	18.102	-17.749	1.00	36.16	C
ATOM	6182	CB GLN A 408	-50.779	19.016	-18.638	1.00	36.63	C
ATOM	6185	CG GLN A 408	-51.625	18.270	-19.672	1.00	38.95	C
ATOM	6188	CD GLN A 408	-52.991	18.930	-19.887	1.00	41.80	C
ATOM	6189	OE1 GLN A 408	-53.739	19.150	-18.921	1.00	43.40	O
ATOM	6190	NE2 GLN A 408	-53.323	19.244	-21.151	1.00	41.68	N
ATOM	6193	C GLN A 408	-48.630	17.725	-18.470	1.00	36.22	C
ATOM	6194	O GLN A 408	-48.456	16.588	-18.920	1.00	35.86	O
ATOM	6196	N ASN A 409	-47.726	18.694	-18.578	1.00	36.64	N
ATOM	6197	CA ASN A 409	-46.400	18.441	-19.136	1.00	36.85	C
ATOM	6199	CB ASN A 409	-46.134	19.306	-20.373	1.00	36.73	C
ATOM	6202	CG ASN A 409	-46.808	18.752	-21.604	1.00	36.19	C
ATOM	6203	OD1 ASN A 409	-46.178	18.072	-22.420	1.00	35.00	O
ATOM	6204	ND2 ASN A 409	-48.110	18.994	-21.720	1.00	34.99	N
ATOM	6207	C ASN A 409	-45.312	18.613	-18.095	1.00	36.96	C
ATOM	6208	O ASN A 409	-45.032	19.720	-17.618	1.00	36.81	O
ATOM	6210	N ILE A 410	-44.713	17.488	-17.734	1.00	37.10	N
ATOM	6211	CA ILE A 410	-43.625	17.508	-16.792	1.00	37.27	C
ATOM	6213	CB ILE A 410	-43.265	16.092	-16.295	1.00	37.38	C
ATOM	6215	CG1 ILE A 410	-42.745	15.209	-17.426	1.00	37.84	C
ATOM	6218	CD1 ILE A 410	-42.416	13.824	-16.985	1.00	38.96	C
ATOM	6222	CG2 ILE A 410	-44.491	15.437	-15.676	1.00	37.67	C
ATOM	6226	C ILE A 410	-42.467	18.179	-17.492	1.00	37.12	C
ATOM	6227	O ILE A 410	-42.150	17.839	-18.617	1.00	36.87	O
ATOM	6229	N LYS A 411	-41.896	19.187	-16.851	1.00	37.38	N
ATOM	6230	CA LYS A 411	-40.680	19.813	-17.332	1.00	37.73	C
ATOM	6232	CB LYS A 411	-40.756	21.332	-17.171	1.00	38.15	C
ATOM	6235	CG LYS A 411	-41.979	21.962	-17.837	1.00	39.89	C
ATOM	6238	CD LYS A 411	-41.797	23.463	-18.080	1.00	42.48	C
ATOM	6241	CE LYS A 411	-43.029	24.069	-18.795	1.00	44.22	C
ATOM	6244	NZ LYS A 411	-42.977	25.576	-18.904	1.00	45.31	N
ATOM	6248	C LYS A 411	-39.562	19.244	-16.489	1.00	37.38	C
ATOM	6249	O LYS A 411	-39.721	19.112	-15.281	1.00	37.28	O
ATOM	6251	N LYS A 412	-38.445	18.886	-17.122	1.00	37.24	N

ATOM	6252	CA	LYS A 412	-37.333	18.243	-16.413	1.00	36.95	C
ATOM	6254	CB	LYS A 412	-36.286	17.701	-17.384	1.00	37.14	C
ATOM	6257	CG	LYS A 412	-35.233	16.818	-16.720	1.00	37.89	C
ATOM	6260	CD	LYS A 412	-34.334	16.125	-17.757	1.00	39.54	C
ATOM	6263	CE	LYS A 412	-33.263	17.068	-18.348	1.00	40.08	C
ATOM	6266	NZ	LYS A 412	-32.161	17.399	-17.376	1.00	40.16	N
ATOM	6270	C	LYS A 412	-36.664	19.177	-15.419	1.00	36.42	C
ATOM	6271	O	LYS A 412	-36.246	18.728	-14.357	1.00	36.56	O
ATOM	6273	N	GLU A 413	-36.570	20.466	-15.744	1.00	35.77	N
ATOM	6274	CA	GLU A 413	-35.948	21.423	-14.820	1.00	35.38	C
ATOM	6276	CB	GLU A 413	-35.688	22.792	-15.489	1.00	35.73	C
ATOM	6279	CG	GLU A 413	-36.510	24.013	-14.973	1.00	37.10	C
ATOM	6282	CD	GLU A 413	-35.621	25.235	-14.633	1.00	38.51	C
ATOM	6283	OE1	GLU A 413	-34.778	25.125	-13.714	1.00	39.24	O
ATOM	6284	OE2	GLU A 413	-35.768	26.306	-15.266	1.00	39.04	O
ATOM	6285	C	GLU A 413	-36.760	21.551	-13.522	1.00	34.43	C
ATOM	6286	O	GLU A 413	-36.196	21.801	-12.460	1.00	34.34	O
ATOM	6288	N	GLU A 414	-38.074	21.360	-13.616	1.00	33.45	N
ATOM	6289	CA	GLU A 414	-38.956	21.407	-12.452	1.00	32.80	C
ATOM	6291	CB	GLU A 414	-40.435	21.425	-12.864	1.00	32.83	C
ATOM	6294	CG	GLU A 414	-40.923	22.736	-13.478	1.00	33.15	C
ATOM	6297	CD	GLU A 414	-42.360	22.668	-14.003	1.00	33.82	C
ATOM	6298	OE1	GLU A 414	-42.876	21.560	-14.231	1.00	35.31	O
ATOM	6299	OE2	GLU A 414	-42.985	23.727	-14.204	1.00	33.98	O
ATOM	6300	C	GLU A 414	-38.715	20.222	-11.536	1.00	32.18	C
ATOM	6301	O	GLU A 414	-38.407	20.409	-10.372	1.00	32.25	O
ATOM	6303	N	ILE A 415	-38.863	19.002	-12.044	1.00	31.63	N
ATOM	6304	CA	ILE A 415	-38.725	17.819	-11.181	1.00	31.28	C
ATOM	6306	CB	ILE A 415	-39.131	16.472	-11.852	1.00	31.14	C
ATOM	6308	CG1	ILE A 415	-38.349	16.204	-13.133	1.00	31.13	C
ATOM	6311	CD1	ILE A 415	-38.786	14.926	-13.836	1.00	30.72	C
ATOM	6315	CG2	ILE A 415	-40.611	16.454	-12.160	1.00	31.12	C
ATOM	6319	C	ILE A 415	-37.316	17.705	-10.636	1.00	31.08	C
ATOM	6320	O	ILE A 415	-37.105	17.188	-9.539	1.00	31.22	O
ATOM	6322	N	GLU A 416	-36.345	18.208	-11.381	1.00	30.67	N
ATOM	6323	CA	GLU A 416	-34.998	18.270	-10.850	1.00	30.55	C
ATOM	6325	CB	GLU A 416	-34.011	18.654	-11.955	1.00	30.91	C
ATOM	6328	CG	GLU A 416	-32.654	17.967	-11.845	1.00	32.55	C
ATOM	6331	CD	GLU A 416	-31.802	18.130	-13.109	1.00	34.91	C
ATOM	6332	OE1	GLU A 416	-32.337	18.563	-14.161	1.00	35.53	O
ATOM	6333	OE2	GLU A 416	-30.590	17.816	-13.052	1.00	36.45	O
ATOM	6334	C	GLU A 416	-34.962	19.249	-9.650	1.00	29.74	C
ATOM	6335	O	GLU A 416	-34.143	19.097	-8.738	1.00	29.73	O
ATOM	6337	N	ASN A 417	-35.864	20.234	-9.652	1.00	28.69	N
ATOM	6338	CA	ASN A 417	-36.060	21.134	-8.503	1.00	27.92	C
ATOM	6340	CB	ASN A 417	-36.632	22.477	-8.963	1.00	27.75	C
ATOM	6343	CG	ASN A 417	-35.572	23.532	-9.095	1.00	26.84	C
ATOM	6344	OD1	ASN A 417	-35.160	24.132	-8.105	1.00	24.88	O
ATOM	6345	ND2	ASN A 417	-35.117	23.765	-10.318	1.00	26.28	N
ATOM	6348	C	ASN A 417	-36.917	20.580	-7.356	1.00	27.41	C
ATOM	6349	O	ASN A 417	-36.650	20.838	-6.187	1.00	27.15	O
ATOM	6351	N	LEU A 418	-37.955	19.834	-7.675	1.00	27.09	N
ATOM	6352	CA	LEU A 418	-38.694	19.156	-6.628	1.00	27.02	C
ATOM	6354	CB	LEU A 418	-39.921	18.423	-7.200	1.00	26.95	C
ATOM	6357	CG	LEU A 418	-41.030	19.284	-7.826	1.00	25.97	C
ATOM	6359	CD1	LEU A 418	-42.127	18.429	-8.447	1.00	24.88	C
ATOM	6363	CD2	LEU A 418	-41.623	20.226	-6.803	1.00	24.50	C
ATOM	6367	C	LEU A 418	-37.763	18.196	-5.860	1.00	27.26	C
ATOM	6368	O	LEU A 418	-37.834	18.123	-4.644	1.00	27.07	O
ATOM	6370	N	GLN A 419	-36.873	17.494	-6.561	1.00	27.67	N

ATOM	6371	CA	GLN A 419	-35.924	16.575	-5.909	1.00	28.28	C
ATOM	6373	CB	GLN A 419	-35.155	15.758	-6.945	1.00	28.36	C
ATOM	6376	CG	GLN A 419	-35.747	14.375	-7.146	1.00	29.43	C
ATOM	6379	CD	GLN A 419	-35.052	13.573	-8.215	1.00	29.76	C
ATOM	6380	OE1	GLN A 419	-34.355	14.112	-9.066	1.00	29.93	O
ATOM	6381	NE2	GLN A 419	-35.252	12.270	-8.184	1.00	31.33	N
ATOM	6384	C	GLN A 419	-34.922	17.219	-4.947	1.00	28.73	C
ATOM	6385	O	GLN A 419	-34.540	16.608	-3.952	1.00	28.90	O
ATOM	6387	N	LYS A 420	-34.490	18.440	-5.240	1.00	29.27	N
ATOM	6388	CA	LYS A 420	-33.608	19.166	-4.335	1.00	29.81	C
ATOM	6390	CB	LYS A 420	-32.758	20.163	-5.129	1.00	30.24	C
ATOM	6393	CG	LYS A 420	-31.592	19.513	-5.867	1.00	31.75	C
ATOM	6396	CD	LYS A 420	-31.259	20.217	-7.183	1.00	34.14	C
ATOM	6399	CE	LYS A 420	-30.140	19.463	-7.953	1.00	36.21	C
ATOM	6402	NZ	LYS A 420	-30.073	19.759	-9.439	1.00	36.76	N
ATOM	6406	C	LYS A 420	-34.375	19.872	-3.204	1.00	29.98	C
ATOM	6407	O	LYS A 420	-33.759	20.547	-2.384	1.00	30.44	O
ATOM	6409	N	TYR A 421	-35.706	19.694	-3.162	1.00	30.07	N
ATOM	6410	CA	TYR A 421	-36.646	20.309	-2.171	1.00	29.73	C
ATOM	6412	CB	TYR A 421	-36.298	19.949	-.712	1.00	29.59	C
ATOM	6415	CG	TYR A 421	-36.386	18.468	-.470	1.00	30.57	C
ATOM	6416	CD1	TYR A 421	-37.582	17.785	-.643	1.00	30.87	C
ATOM	6418	CE1	TYR A 421	-37.665	16.414	-.448	1.00	31.28	C
ATOM	6420	CZ	TYR A 421	-36.549	15.705	-.068	1.00	31.90	C
ATOM	6421	OH	TYR A 421	-36.620	14.347	.124	1.00	32.43	O
ATOM	6423	CE2	TYR A 421	-35.355	16.354	.119	1.00	32.68	C
ATOM	6425	CD2	TYR A 421	-35.274	17.736	-.095	1.00	32.42	C
ATOM	6427	C	TYR A 421	-36.850	21.816	-2.347	1.00	29.24	C
ATOM	6428	O	TYR A 421	-36.709	22.590	-1.401	1.00	29.26	O
ATOM	6430	N	HIS A 422	-37.219	22.214	-3.563	1.00	28.65	N
ATOM	6431	CA	HIS A 422	-37.540	23.604	-3.859	1.00	28.26	C
ATOM	6433	CB	HIS A 422	-38.150	23.715	-5.256	1.00	28.33	C
ATOM	6436	CG	HIS A 422	-38.299	25.124	-5.748	1.00	28.44	C
ATOM	6437	ND1	HIS A 422	-37.225	25.904	-6.106	1.00	28.97	N
ATOM	6439	CE1	HIS A 422	-37.660	27.082	-6.520	1.00	29.40	C
ATOM	6441	NE2	HIS A 422	-38.977	27.090	-6.443	1.00	28.14	N
ATOM	6443	CD2	HIS A 422	-39.399	25.879	-5.965	1.00	27.76	C
ATOM	6445	C	HIS A 422	-38.511	24.147	-2.816	1.00	27.74	C
ATOM	6446	O	HIS A 422	-39.357	23.408	-2.310	1.00	27.53	O
ATOM	6448	N	ASP A 423	-38.384	25.437	-2.505	1.00	27.12	N
ATOM	6449	CA	ASP A 423	-39.186	26.072	-1.454	1.00	26.63	C
ATOM	6451	CB	ASP A 423	-38.877	27.575	-1.349	1.00	27.05	C
ATOM	6454	CG	ASP A 423	-37.456	27.875	-.851	1.00	28.35	C
ATOM	6455	OD1	ASP A 423	-36.951	27.121	.018	1.00	30.99	O
ATOM	6456	OD2	ASP A 423	-36.860	28.884	-1.322	1.00	27.89	O
ATOM	6457	C	ASP A 423	-40.689	25.883	-1.698	1.00	25.71	C
ATOM	6458	O	ASP A 423	-41.487	25.831	-.749	1.00	26.26	O
ATOM	6460	N	ILE A 424	-41.075	25.781	-2.966	1.00	24.01	N
ATOM	6461	CA	ILE A 424	-42.479	25.596	-3.327	1.00	22.72	C
ATOM	6463	CB	ILE A 424	-42.639	25.391	-4.843	1.00	22.40	C
ATOM	6465	CG1	ILE A 424	-44.029	25.668	-5.318	1.00	21.78	C
ATOM	6468	CD1	ILE A 424	-44.086	25.499	-6.793	1.00	22.22	C
ATOM	6472	CG2	ILE A 424	-42.340	23.985	-5.255	1.00	22.87	C
ATOM	6476	C	ILE A 424	-43.103	24.442	-2.562	1.00	21.91	C
ATOM	6477	O	ILE A 424	-44.238	24.574	-2.115	1.00	21.73	O
ATOM	6479	N	ILE A 425	-42.365	23.334	-2.399	1.00	21.23	N
ATOM	6480	CA	ILE A 425	-42.853	22.173	-1.634	1.00	20.73	C
ATOM	6482	CB	ILE A 425	-42.622	20.820	-2.330	1.00	20.10	C
ATOM	6484	CG1	ILE A 425	-41.158	20.409	-2.290	1.00	18.72	C
ATOM	6487	CD1	ILE A 425	-40.923	19.075	-2.948	1.00	18.27	C

ATOM	6491	CG2 ILE A 425	-43.146	20.840	-3.741	1.00	19.57	C
ATOM	6495	C ILE A 425	-42.252	22.060	-.242	1.00	21.21	C
ATOM	6496	O ILE A 425	-42.810	21.370	.613	1.00	21.30	O
ATOM	6498	N SER A 426	-41.129	22.721	-.002	1.00	21.56	N
ATOM	6499	CA SER A 426	-40.492	22.620	1.299	1.00	22.23	C
ATOM	6501	CB SER A 426	-39.043	23.082	1.234	1.00	22.36	C
ATOM	6504	OG SER A 426	-38.977	24.501	1.342	1.00	24.09	O
ATOM	6506	C SER A 426	-41.242	23.460	2.334	1.00	22.39	C
ATOM	6507	O SER A 426	-41.390	23.051	3.491	1.00	22.95	O
ATOM	6509	N ARG A 427	-41.701	24.643	1.937	1.00	22.08	N
ATOM	6510	CA ARG A 427	-42.332	25.541	2.905	1.00	21.87	C
ATOM	6512	CB ARG A 427	-42.466	26.955	2.345	1.00	21.99	C
ATOM	6515	CG ARG A 427	-41.170	27.696	2.493	1.00	23.36	C
ATOM	6518	CD ARG A 427	-40.954	28.727	1.433	1.00	26.07	C
ATOM	6521	NE ARG A 427	-39.680	29.407	1.667	1.00	28.18	N
ATOM	6523	CZ ARG A 427	-39.127	30.298	.844	1.00	29.21	C
ATOM	6524	NH1 ARG A 427	-39.723	30.633	-.298	1.00	28.60	N
ATOM	6527	NH2 ARG A 427	-37.962	30.854	1.173	1.00	30.48	N
ATOM	6530	C ARG A 427	-43.651	25.003	3.434	1.00	21.21	C
ATOM	6531	O ARG A 427	-43.817	24.891	4.645	1.00	21.11	O
ATOM	6533	N PRO A 428	-44.573	24.622	2.541	1.00	20.52	N
ATOM	6534	CA PRO A 428	-45.789	23.999	3.040	1.00	20.27	C
ATOM	6536	CB PRO A 428	-46.410	23.369	1.791	1.00	20.44	C
ATOM	6539	CG PRO A 428	-45.864	24.119	.666	1.00	20.63	C
ATOM	6542	CD PRO A 428	-44.521	24.644	1.074	1.00	20.54	C
ATOM	6545	C PRO A 428	-45.450	22.917	4.050	1.00	19.75	C
ATOM	6546	O PRO A 428	-46.089	22.834	5.095	1.00	19.77	O
ATOM	6547	N SER A 429	-44.424	22.123	3.747	1.00	18.97	N
ATOM	6548	CA SER A 429	-44.008	21.045	4.636	1.00	18.45	C
ATOM	6550	CB SER A 429	-42.954	20.180	3.972	1.00	18.35	C
ATOM	6553	OG SER A 429	-43.452	19.716	2.733	1.00	19.12	O
ATOM	6555	C SER A 429	-43.519	21.525	5.985	1.00	17.97	C
ATOM	6556	O SER A 429	-43.733	20.846	6.968	1.00	17.88	O
ATOM	6558	N HIS A 430	-42.879	22.689	6.053	1.00	17.73	N
ATOM	6559	CA HIS A 430	-42.593	23.289	7.359	1.00	17.71	C
ATOM	6561	CB HIS A 430	-41.937	24.668	7.241	1.00	17.97	C
ATOM	6564	CG HIS A 430	-40.558	24.663	6.651	1.00	18.87	C
ATOM	6565	ND1 HIS A 430	-39.940	23.524	6.186	1.00	19.85	N
ATOM	6567	CE1 HIS A 430	-38.746	23.836	5.709	1.00	19.43	C
ATOM	6569	NE2 HIS A 430	-38.573	25.137	5.838	1.00	19.06	N
ATOM	6571	CD2 HIS A 430	-39.693	25.681	6.418	1.00	19.32	C
ATOM	6573	C HIS A 430	-43.914	23.444	8.131	1.00	17.31	C
ATOM	6574	O HIS A 430	-44.023	23.013	9.287	1.00	17.06	O
ATOM	6576	N ILE A 431	-44.913	24.050	7.475	1.00	16.74	N
ATOM	6577	CA ILE A 431	-46.218	24.304	8.091	1.00	16.15	C
ATOM	6579	CB ILE A 431	-47.174	25.057	7.174	1.00	16.18	C
ATOM	6581	CG1 ILE A 431	-46.613	26.428	6.801	1.00	17.05	C
ATOM	6584	CD1 ILE A 431	-46.441	27.352	7.994	1.00	18.34	C
ATOM	6588	CG2 ILE A 431	-48.496	25.263	7.860	1.00	15.28	C
ATOM	6592	C ILE A 431	-46.899	23.021	8.476	1.00	15.75	C
ATOM	6593	O ILE A 431	-47.624	22.980	9.445	1.00	15.95	O
ATOM	6595	N PHE A 432	-46.663	21.969	7.714	1.00	15.59	N
ATOM	6596	CA PHE A 432	-47.192	20.652	8.041	1.00	15.52	C
ATOM	6598	CB PHE A 432	-46.837	19.683	6.913	1.00	15.23	C
ATOM	6601	CG PHE A 432	-47.451	18.318	7.031	1.00	15.98	C
ATOM	6602	CD1 PHE A 432	-48.324	17.972	8.052	1.00	16.55	C
ATOM	6604	CE1 PHE A 432	-48.867	16.700	8.102	1.00	16.04	C
ATOM	6606	CZ PHE A 432	-48.559	15.774	7.130	1.00	15.29	C
ATOM	6608	CE2 PHE A 432	-47.718	16.107	6.114	1.00	15.37	C
ATOM	6610	CD2 PHE A 432	-47.172	17.364	6.061	1.00	16.74	C

ATOM	6612	C	PHE A 432	-46.640	20.191	9.392	1.00	15.54	C
ATOM	6613	O	PHE A 432	-47.383	20.006	10.343	1.00	15.37	O
ATOM	6615	N	ARG A 433	-45.328	20.041	9.483	1.00	15.98	N
ATOM	6616	CA	ARG A 433	-44.695	19.542	10.709	1.00	16.20	C
ATOM	6618	CB	ARG A 433	-43.176	19.476	10.526	1.00	15.85	C
ATOM	6621	CG	ARG A 433	-42.411	19.183	11.804	1.00	16.11	C
ATOM	6624	CD	ARG A 433	-42.848	17.890	12.489	1.00	15.70	C
ATOM	6627	NE	ARG A 433	-42.811	16.747	11.588	1.00	15.64	N
ATOM	6629	CZ	ARG A 433	-43.377	15.571	11.838	1.00	16.27	C
ATOM	6630	NH1	ARG A 433	-44.029	15.359	12.965	1.00	16.82	N
ATOM	6633	NH2	ARG A 433	-43.302	14.598	10.948	1.00	16.96	N
ATOM	6636	C	ARG A 433	-45.029	20.411	11.935	1.00	16.46	C
ATOM	6637	O	ARG A 433	-45.296	19.890	13.031	1.00	16.31	O
ATOM	6639	N	LEU A 434	-45.006	21.731	11.724	1.00	16.47	N
ATOM	6640	CA	LEU A 434	-45.123	22.701	12.800	1.00	15.98	C
ATOM	6642	CB	LEU A 434	-44.771	24.114	12.303	1.00	15.83	C
ATOM	6645	CG	LEU A 434	-43.287	24.486	12.218	1.00	14.95	C
ATOM	6647	CD1	LEU A 434	-43.130	25.918	11.826	1.00	15.09	C
ATOM	6651	CD2	LEU A 434	-42.600	24.281	13.536	1.00	14.37	C
ATOM	6655	C	LEU A 434	-46.517	22.660	13.387	1.00	16.19	C
ATOM	6656	O	LEU A 434	-46.669	22.564	14.604	1.00	15.92	O
ATOM	6658	N	CYS A 435	-47.529	22.715	12.520	1.00	16.78	N
ATOM	6659	CA	CYS A 435	-48.936	22.591	12.940	1.00	17.48	C
ATOM	6661	CB	CYS A 435	-49.870	22.617	11.726	1.00	17.54	C
ATOM	6664	SG	CYS A 435	-50.110	24.233	10.985	1.00	18.80	S
ATOM	6666	C	CYS A 435	-49.192	21.296	13.708	1.00	17.84	C
ATOM	6667	O	CYS A 435	-50.043	21.250	14.597	1.00	17.92	O
ATOM	6669	N	ASN A 436	-48.446	20.257	13.333	1.00	18.24	N
ATOM	6670	CA	ASN A 436	-48.603	18.913	13.846	1.00	18.50	C
ATOM	6672	CB	ASN A 436	-47.936	17.949	12.867	1.00	18.63	C
ATOM	6675	CG	ASN A 436	-48.156	16.487	13.213	1.00	18.58	C
ATOM	6676	OD1	ASN A 436	-48.810	16.129	14.185	1.00	18.80	O
ATOM	6677	ND2	ASN A 436	-47.588	15.633	12.398	1.00	19.29	N
ATOM	6680	C	ASN A 436	-47.975	18.755	15.210	1.00	18.83	C
ATOM	6681	O	ASN A 436	-48.551	18.153	16.115	1.00	19.00	O
ATOM	6683	N	ASP A 437	-46.763	19.255	15.351	1.00	19.29	N
ATOM	6684	CA	ASP A 437	-46.097	19.179	16.634	1.00	19.62	C
ATOM	6686	CB	ASP A 437	-44.615	19.522	16.510	1.00	19.45	C
ATOM	6689	CG	ASP A 437	-43.849	18.480	15.734	1.00	19.36	C
ATOM	6690	OD1	ASP A 437	-44.491	17.586	15.158	1.00	19.89	O
ATOM	6691	OD2	ASP A 437	-42.608	18.543	15.699	1.00	19.86	O
ATOM	6692	C	ASP A 437	-46.799	20.106	17.597	1.00	20.01	C
ATOM	6693	O	ASP A 437	-46.916	19.790	18.779	1.00	20.30	O
ATOM	6695	N	LEU A 438	-47.288	21.233	17.092	1.00	20.40	N
ATOM	6696	CA	LEU A 438	-47.996	22.174	17.940	1.00	21.03	C
ATOM	6698	CB	LEU A 438	-48.510	23.363	17.126	1.00	21.01	C
ATOM	6701	CG	LEU A 438	-47.575	24.565	17.078	1.00	20.32	C
ATOM	6703	CD1	LEU A 438	-47.938	25.524	15.952	1.00	19.14	C
ATOM	6707	CD2	LEU A 438	-47.611	25.263	18.415	1.00	19.73	C
ATOM	6711	C	LEU A 438	-49.153	21.496	18.679	1.00	21.81	C
ATOM	6712	O	LEU A 438	-49.277	21.633	19.902	1.00	21.55	O
ATOM	6714	N	ALA A 439	-49.973	20.754	17.932	1.00	22.76	N
ATOM	6715	CA	ALA A 439	-51.138	20.055	18.486	1.00	23.58	C
ATOM	6717	CB	ALA A 439	-51.905	19.360	17.386	1.00	23.57	C
ATOM	6721	C	ALA A 439	-50.746	19.042	19.542	1.00	24.47	C
ATOM	6722	O	ALA A 439	-51.386	18.951	20.575	1.00	24.63	O
ATOM	6724	N	SER A 440	-49.685	18.291	19.276	1.00	25.78	N
ATOM	6725	CA	SER A 440	-49.242	17.227	20.168	1.00	26.84	C
ATOM	6727	CB	SER A 440	-48.566	16.130	19.353	1.00	26.92	C
ATOM	6730	OG	SER A 440	-47.321	16.582	18.859	1.00	27.46	O

ATOM	6732	C	SER A 440	-48.270	17.691	21.244	1.00	27.68	C
ATOM	6733	O	SER A 440	-47.934	16.923	22.132	1.00	28.00	O
ATOM	6735	N	ALA A 441	-47.822	18.935	21.174	1.00	28.83	N
ATOM	6736	CA	ALA A 441	-46.760	19.419	22.056	1.00	29.73	C
ATOM	6738	CB	ALA A 441	-46.529	20.910	21.841	1.00	29.84	C
ATOM	6742	C	ALA A 441	-46.982	19.132	23.539	1.00	30.57	C
ATOM	6743	O	ALA A 441	-46.306	18.274	24.093	1.00	30.50	O
ATOM	6745	N	SER A 442	-47.929	19.827	24.172	1.00	31.90	N
ATOM	6746	CA	SER A 442	-47.999	19.862	25.654	1.00	32.95	C
ATOM	6748	CB	SER A 442	-49.043	20.870	26.169	1.00	32.91	C
ATOM	6751	OG	SER A 442	-50.331	20.614	25.650	1.00	33.50	O
ATOM	6753	C	SER A 442	-48.219	18.497	26.289	1.00	33.69	C
ATOM	6754	O	SER A 442	-47.754	18.246	27.397	1.00	33.59	O
ATOM	6756	N	ALA A 443	-48.922	17.623	25.578	1.00	34.93	N
ATOM	6757	CA	ALA A 443	-49.050	16.231	25.980	1.00	35.81	C
ATOM	6759	CB	ALA A 443	-49.981	15.488	25.030	1.00	35.72	C
ATOM	6763	C	ALA A 443	-47.663	15.580	25.996	1.00	36.71	C
ATOM	6764	O	ALA A 443	-47.198	15.115	27.047	1.00	36.94	O
ATOM	6766	N	GLU A 444	-46.999	15.567	24.838	1.00	37.47	N
ATOM	6767	CA	GLU A 444	-45.688	14.930	24.718	1.00	37.93	C
ATOM	6769	CB	GLU A 444	-45.164	14.996	23.277	1.00	37.99	C
ATOM	6772	CG	GLU A 444	-45.952	14.100	22.326	1.00	39.24	C
ATOM	6775	CD	GLU A 444	-45.419	14.080	20.886	1.00	41.03	C
ATOM	6776	OE1	GLU A 444	-44.612	14.958	20.510	1.00	42.16	O
ATOM	6777	OE2	GLU A 444	-45.827	13.178	20.119	1.00	42.36	O
ATOM	6778	C	GLU A 444	-44.694	15.534	25.708	1.00	38.09	C
ATOM	6779	O	GLU A 444	-43.924	14.803	26.324	1.00	38.17	O
ATOM	6781	N	ILE A 445	-44.739	16.851	25.889	1.00	38.41	N
ATOM	6782	CA	ILE A 445	-43.829	17.530	26.814	1.00	38.75	C
ATOM	6784	CB	ILE A 445	-43.802	19.065	26.587	1.00	38.69	C
ATOM	6786	CG1	ILE A 445	-43.314	19.396	25.170	1.00	38.18	C
ATOM	6789	CD1	ILE A 445	-43.702	20.775	24.712	1.00	37.58	C
ATOM	6793	CG2	ILE A 445	-42.909	19.753	27.617	1.00	38.26	C
ATOM	6797	C	ILE A 445	-44.243	17.194	28.246	1.00	39.35	C
ATOM	6798	O	ILE A 445	-45.039	17.901	28.870	1.00	39.44	O
ATOM	6800	N	ALA A 446	-43.696	16.092	28.750	1.00	40.02	N
ATOM	6801	CA	ALA A 446	-44.082	15.540	30.045	1.00	40.51	C
ATOM	6803	CB	ALA A 446	-43.815	16.552	31.170	1.00	40.55	C
ATOM	6807	C	ALA A 446	-45.559	15.130	30.029	1.00	40.86	C
ATOM	6808	O	ALA A 446	-46.407	16.021	30.121	1.00	40.77	O
ATOM	6810	N	ARG A 447	-45.925	13.842	29.899	1.00	41.31	N
ATOM	6811	CA	ARG A 447	-45.072	12.633	29.691	1.00	41.58	C
ATOM	6813	CB	ARG A 447	-45.272	12.092	28.258	1.00	41.84	C
ATOM	6816	CG	ARG A 447	-46.667	11.495	27.991	1.00	42.69	C
ATOM	6819	CD	ARG A 447	-46.778	10.901	26.582	1.00	43.92	C
ATOM	6822	NE	ARG A 447	-47.764	11.589	25.744	1.00	45.34	N
ATOM	6824	CZ	ARG A 447	-47.887	11.435	24.420	1.00	46.56	C
ATOM	6825	NH1	ARG A 447	-47.077	10.623	23.742	1.00	46.71	N
ATOM	6828	NH2	ARG A 447	-48.828	12.107	23.757	1.00	46.97	N
ATOM	6831	C	ARG A 447	-43.581	12.728	30.054	1.00	41.38	C
ATOM	6832	O	ARG A 447	-43.231	12.747	31.238	1.00	41.92	O
ATOM	6834	N	GLY A 448	-42.710	12.736	29.051	1.00	40.88	N
ATOM	6835	CA	GLY A 448	-41.312	13.128	29.234	1.00	40.44	C
ATOM	6838	C	GLY A 448	-40.579	13.316	27.915	1.00	40.08	C
ATOM	6839	O	GLY A 448	-39.364	13.488	27.901	1.00	40.05	O
ATOM	6841	N	GLU A 449	-41.323	13.332	26.809	1.00	39.56	N
ATOM	6842	CA	GLU A 449	-40.750	13.132	25.486	1.00	39.19	C
ATOM	6844	CB	GLU A 449	-41.827	12.668	24.502	1.00	39.52	C
ATOM	6847	CG	GLU A 449	-42.403	11.290	24.810	1.00	40.61	C
ATOM	6850	CD	GLU A 449	-43.348	10.782	23.715	1.00	42.21	C

ATOM	6851	OE1	GLU A 449	-43.068	11.019	22.503	1.00	41.83	O
ATOM	6852	OE2	GLU A 449	-44.368	10.141	24.082	1.00	42.41	O
ATOM	6853	C	GLU A 449	-40.037	14.365	24.921	1.00	38.34	C
ATOM	6854	O	GLU A 449	-40.511	15.497	25.031	1.00	38.01	O
ATOM	6856	N	THR A 450	-38.904	14.106	24.279	1.00	37.37	N
ATOM	6857	CA	THR A 450	-38.047	15.138	23.726	1.00	36.37	C
ATOM	6859	CB	THR A 450	-36.607	14.985	24.285	1.00	36.39	C
ATOM	6861	OG1	THR A 450	-35.846	16.153	23.982	1.00	36.63	O
ATOM	6863	CG2	THR A 450	-35.888	13.749	23.706	1.00	36.52	C
ATOM	6867	C	THR A 450	-38.013	15.100	22.190	1.00	35.36	C
ATOM	6868	O	THR A 450	-37.175	15.764	21.583	1.00	35.43	O
ATOM	6870	N	ALA A 451	-38.917	14.338	21.563	1.00	33.98	N
ATOM	6871	CA	ALA A 451	-38.920	14.182	20.096	1.00	32.81	C
ATOM	6873	CB	ALA A 451	-39.030	12.708	19.708	1.00	32.99	C
ATOM	6877	C	ALA A 451	-40.043	14.987	19.459	1.00	31.46	C
ATOM	6878	O	ALA A 451	-41.050	14.432	18.996	1.00	31.06	O
ATOM	6880	N	ASN A 452	-39.848	16.302	19.432	1.00	29.88	N
ATOM	6881	CA	ASN A 452	-40.881	17.225	18.973	1.00	28.81	C
ATOM	6883	CB	ASN A 452	-41.957	17.329	20.051	1.00	28.57	C
ATOM	6886	CG	ASN A 452	-42.969	18.380	19.757	1.00	28.30	C
ATOM	6887	OD1	ASN A 452	-42.622	19.517	19.468	1.00	29.35	O
ATOM	6888	ND2	ASN A 452	-44.236	18.018	19.836	1.00	28.36	N
ATOM	6891	C	ASN A 452	-40.285	18.590	18.641	1.00	27.76	C
ATOM	6892	O	ASN A 452	-39.365	19.018	19.296	1.00	28.04	O
ATOM	6894	N	SER A 453	-40.807	19.269	17.630	1.00	26.77	N
ATOM	6895	CA	SER A 453	-40.226	20.533	17.182	1.00	26.38	C
ATOM	6897	CB	SER A 453	-40.912	21.015	15.903	1.00	26.40	C
ATOM	6900	OG	SER A 453	-40.796	20.044	14.876	1.00	25.79	O
ATOM	6902	C	SER A 453	-40.253	21.643	18.233	1.00	26.12	C
ATOM	6903	O	SER A 453	-39.280	22.361	18.385	1.00	25.97	O
ATOM	6905	N	VAL A 454	-41.368	21.784	18.942	1.00	26.15	N
ATOM	6906	CA	VAL A 454	-41.518	22.780	20.022	1.00	26.21	C
ATOM	6908	CB	VAL A 454	-42.975	22.818	20.532	1.00	26.04	C
ATOM	6910	CG1	VAL A 454	-43.122	23.787	21.694	1.00	25.24	C
ATOM	6914	CG2	VAL A 454	-43.913	23.190	19.401	1.00	26.32	C
ATOM	6918	C	VAL A 454	-40.617	22.451	21.211	1.00	26.63	C
ATOM	6919	O	VAL A 454	-40.173	23.328	21.959	1.00	26.53	O
ATOM	6921	N	SER A 455	-40.374	21.160	21.380	1.00	27.21	N
ATOM	6922	CA	SER A 455	-39.541	20.651	22.437	1.00	27.60	C
ATOM	6924	CB	SER A 455	-39.677	19.141	22.490	1.00	27.31	C
ATOM	6927	OG	SER A 455	-38.922	18.625	23.545	1.00	28.13	O
ATOM	6929	C	SER A 455	-38.096	21.045	22.181	1.00	28.29	C
ATOM	6930	O	SER A 455	-37.445	21.628	23.043	1.00	28.63	O
ATOM	6932	N	CYS A 456	-37.599	20.748	20.988	1.00	29.04	N
ATOM	6933	CA	CYS A 456	-36.219	21.082	20.640	1.00	29.65	C
ATOM	6935	CB	CYS A 456	-35.868	20.535	19.256	1.00	29.52	C
ATOM	6938	SG	CYS A 456	-35.959	18.734	19.163	1.00	29.71	S
ATOM	6940	C	CYS A 456	-35.955	22.589	20.710	1.00	30.26	C
ATOM	6941	O	CYS A 456	-34.887	23.007	21.162	1.00	30.49	O
ATOM	6943	N	TYR A 457	-36.924	23.399	20.283	1.00	30.88	N
ATOM	6944	CA	TYR A 457	-36.767	24.850	20.308	1.00	31.45	C
ATOM	6946	CB	TYR A 457	-37.967	25.557	19.648	1.00	31.55	C
ATOM	6949	CG	TYR A 457	-37.691	26.991	19.178	1.00	32.09	C
ATOM	6950	CD1	TYR A 457	-37.378	27.269	17.843	1.00	31.69	C
ATOM	6952	CE1	TYR A 457	-37.128	28.569	17.418	1.00	31.73	C
ATOM	6954	CZ	TYR A 457	-37.186	29.613	18.328	1.00	32.84	C
ATOM	6955	OH	TYR A 457	-36.945	30.917	17.927	1.00	33.81	O
ATOM	6957	CE2	TYR A 457	-37.492	29.364	19.658	1.00	32.72	C
ATOM	6959	CD2	TYR A 457	-37.746	28.064	20.074	1.00	32.44	C
ATOM	6961	C	TYR A 457	-36.571	25.295	21.758	1.00	31.99	C

ATOM	6962	O	TYR A 457	-35.648	26.053	22.052	1.00	31.98	O
ATOM	6964	N	MET A 458	-37.419	24.799	22.661	1.00	32.76	N
ATOM	6965	CA	MET A 458	-37.211	24.987	24.105	1.00	33.37	C
ATOM	6967	CB	MET A 458	-38.151	24.097	24.921	1.00	33.43	C
ATOM	6970	CG	MET A 458	-39.570	24.593	25.059	1.00	33.70	C
ATOM	6973	SD	MET A 458	-40.641	23.312	25.737	1.00	34.59	S
ATOM	6974	CE	MET A 458	-39.709	22.763	27.179	1.00	34.71	C
ATOM	6978	C	MET A 458	-35.784	24.641	24.517	1.00	33.78	C
ATOM	6979	O	MET A 458	-35.101	25.437	25.151	1.00	33.73	O
ATOM	6981	N	ARG A 459	-35.347	23.445	24.145	1.00	34.46	N
ATOM	6982	CA	ARG A 459	-34.051	22.926	24.581	1.00	35.15	C
ATOM	6984	CB	ARG A 459	-33.876	21.447	24.171	1.00	35.68	C
ATOM	6987	CG	ARG A 459	-32.494	20.864	24.500	1.00	37.23	C
ATOM	6990	CD	ARG A 459	-32.491	19.345	24.752	1.00	39.80	C
ATOM	6993	NE	ARG A 459	-33.140	18.538	23.711	1.00	42.88	N
ATOM	6995	CZ	ARG A 459	-32.704	18.407	22.453	1.00	45.43	C
ATOM	6996	NH1	ARG A 459	-31.615	19.062	22.035	1.00	47.50	N
ATOM	6999	NH2	ARG A 459	-33.366	17.629	21.595	1.00	44.77	N
ATOM	7002	C	ARG A 459	-32.878	23.749	24.072	1.00	34.92	C
ATOM	7003	O	ARG A 459	-32.007	24.109	24.849	1.00	35.06	O
ATOM	7005	N	THR A 460	-32.854	24.047	22.778	1.00	34.96	N
ATOM	7006	CA	THR A 460	-31.674	24.678	22.163	1.00	34.91	C
ATOM	7008	CB	THR A 460	-31.494	24.263	20.680	1.00	34.88	C
ATOM	7010	OG1	THR A 460	-32.330	25.067	19.841	1.00	34.39	O
ATOM	7012	CG2	THR A 460	-31.825	22.776	20.488	1.00	35.44	C
ATOM	7016	C	THR A 460	-31.672	26.210	22.258	1.00	34.78	C
ATOM	7017	O	THR A 460	-30.673	26.837	21.911	1.00	34.92	O
ATOM	7019	N	LYS A 461	-32.781	26.806	22.701	1.00	34.41	N
ATOM	7020	CA	LYS A 461	-32.811	28.231	23.039	1.00	34.25	C
ATOM	7022	CB	LYS A 461	-33.921	28.947	22.265	1.00	34.56	C
ATOM	7025	CG	LYS A 461	-33.717	28.986	20.750	1.00	35.84	C
ATOM	7028	CD	LYS A 461	-32.519	29.859	20.347	1.00	37.64	C
ATOM	7031	CE	LYS A 461	-32.232	29.804	18.838	1.00	38.44	C
ATOM	7034	NZ	LYS A 461	-33.141	30.671	18.027	1.00	38.31	N
ATOM	7038	C	LYS A 461	-32.974	28.476	24.544	1.00	33.73	C
ATOM	7039	O	LYS A 461	-32.994	29.626	24.983	1.00	33.23	O
ATOM	7041	N	GLY A 462	-33.078	27.395	25.321	1.00	33.48	N
ATOM	7042	CA	GLY A 462	-33.173	27.464	26.784	1.00	33.23	C
ATOM	7045	C	GLY A 462	-34.332	28.320	27.239	1.00	32.95	C
ATOM	7046	O	GLY A 462	-34.144	29.258	28.004	1.00	33.17	O
ATOM	7048	N	ILE A 463	-35.528	28.000	26.757	1.00	32.52	N
ATOM	7049	CA	ILE A 463	-36.694	28.843	26.971	1.00	32.32	C
ATOM	7051	CB	ILE A 463	-37.006	29.724	25.734	1.00	32.47	C
ATOM	7053	CG1	ILE A 463	-37.152	28.866	24.463	1.00	32.41	C
ATOM	7056	CD1	ILE A 463	-37.359	29.675	23.191	1.00	32.02	C
ATOM	7060	CG2	ILE A 463	-35.937	30.810	25.571	1.00	32.44	C
ATOM	7064	C	ILE A 463	-37.902	28.006	27.309	1.00	32.15	C
ATOM	7065	O	ILE A 463	-37.886	26.799	27.145	1.00	32.01	O
ATOM	7067	N	SER A 464	-38.950	28.670	27.780	1.00	32.27	N
ATOM	7068	CA	SER A 464	-40.163	28.006	28.239	1.00	32.52	C
ATOM	7070	CB	SER A 464	-40.964	28.974	29.110	1.00	32.65	C
ATOM	7073	OG	SER A 464	-41.112	30.224	28.457	1.00	32.69	O
ATOM	7075	C	SER A 464	-41.036	27.522	27.079	1.00	32.53	C
ATOM	7076	O	SER A 464	-40.968	28.069	25.986	1.00	32.62	O
ATOM	7078	N	GLU A 465	-41.857	26.500	27.332	1.00	32.47	N
ATOM	7079	CA	GLU A 465	-42.836	26.010	26.353	1.00	32.33	C
ATOM	7081	CB	GLU A 465	-43.689	24.868	26.935	1.00	32.34	C
ATOM	7084	CG	GLU A 465	-44.756	24.312	25.974	1.00	32.37	C
ATOM	7087	CD	GLU A 465	-45.675	23.279	26.611	1.00	32.47	C
ATOM	7088	OE1	GLU A 465	-45.297	22.659	27.630	1.00	31.63	O

ATOM	7089	OE2	GLU A 465	-46.785	23.086	26.074	1.00	32.77	O
ATOM	7090	C	GLU A 465	-43.758	27.124	25.864	1.00	32.33	C
ATOM	7091	O	GLU A 465	-44.078	27.183	24.676	1.00	32.60	O
ATOM	7093	N	GLU A 466	-44.196	28.001	26.767	1.00	32.08	N
ATOM	7094	CA	GLU A 466	-45.085	29.097	26.380	1.00	31.79	C
ATOM	7096	CB	GLU A 466	-45.606	29.832	27.624	1.00	31.91	C
ATOM	7099	CG	GLU A 466	-46.503	31.041	27.304	1.00	32.82	C
ATOM	7102	CD	GLU A 466	-47.426	31.472	28.452	1.00	33.67	C
ATOM	7103	OE1	GLU A 466	-47.519	30.756	29.478	1.00	35.09	O
ATOM	7104	OE2	GLU A 466	-48.077	32.533	28.312	1.00	33.12	O
ATOM	7105	C	GLU A 466	-44.392	30.052	25.386	1.00	31.13	C
ATOM	7106	O	GLU A 466	-45.032	30.603	24.492	1.00	30.87	O
ATOM	7108	N	LEU A 467	-43.077	30.194	25.531	1.00	30.57	N
ATOM	7109	CA	LEU A 467	-42.265	31.102	24.708	1.00	30.18	C
ATOM	7111	CB	LEU A 467	-41.055	31.604	25.524	1.00	30.38	C
ATOM	7114	CG	LEU A 467	-40.674	33.087	25.463	1.00	30.39	C
ATOM	7116	CD1	LEU A 467	-41.608	33.900	26.362	1.00	30.23	C
ATOM	7120	CD2	LEU A 467	-39.214	33.308	25.872	1.00	30.31	C
ATOM	7124	C	LEU A 467	-41.773	30.401	23.439	1.00	29.47	C
ATOM	7125	O	LEU A 467	-41.660	31.010	22.382	1.00	29.17	O
ATOM	7127	N	ALA A 468	-41.445	29.121	23.562	1.00	28.93	N
ATOM	7128	CA	ALA A 468	-41.095	28.310	22.404	1.00	28.45	C
ATOM	7130	CB	ALA A 468	-40.613	26.935	22.830	1.00	28.10	C
ATOM	7134	C	ALA A 468	-42.306	28.193	21.494	1.00	28.08	C
ATOM	7135	O	ALA A 468	-42.179	28.321	20.285	1.00	28.24	O
ATOM	7137	N	THR A 469	-43.479	27.964	22.080	1.00	27.58	N
ATOM	7138	CA	THR A 469	-44.719	27.870	21.322	1.00	27.26	C
ATOM	7140	CB	THR A 469	-45.928	27.710	22.261	1.00	27.28	C
ATOM	7142	OG1	THR A 469	-46.063	26.330	22.627	1.00	27.55	O
ATOM	7144	CG2	THR A 469	-47.222	28.185	21.597	1.00	27.34	C
ATOM	7148	C	THR A 469	-44.916	29.099	20.449	1.00	27.10	C
ATOM	7149	O	THR A 469	-45.176	28.989	19.250	1.00	27.11	O
ATOM	7151	N	GLU A 470	-44.774	30.267	21.061	1.00	26.88	N
ATOM	7152	CA	GLU A 470	-44.944	31.542	20.376	1.00	26.84	C
ATOM	7154	CB	GLU A 470	-44.832	32.663	21.400	1.00	27.07	C
ATOM	7157	CG	GLU A 470	-45.357	34.003	20.958	1.00	28.47	C
ATOM	7160	CD	GLU A 470	-45.140	35.087	22.008	1.00	30.51	C
ATOM	7161	OE1	GLU A 470	-45.029	34.775	23.221	1.00	31.53	O
ATOM	7162	OE2	GLU A 470	-45.083	36.269	21.608	1.00	32.59	O
ATOM	7163	C	GLU A 470	-43.926	31.750	19.239	1.00	26.48	C
ATOM	7164	O	GLU A 470	-44.261	32.296	18.184	1.00	26.29	O
ATOM	7166	N	SER A 471	-42.686	31.316	19.452	1.00	26.12	N
ATOM	7167	CA	SER A 471	-41.660	31.386	18.408	1.00	25.85	C
ATOM	7169	CB	SER A 471	-40.296	30.933	18.941	1.00	25.96	C
ATOM	7172	OG	SER A 471	-39.656	31.939	19.711	1.00	26.53	O
ATOM	7174	C	SER A 471	-42.036	30.529	17.208	1.00	25.35	C
ATOM	7175	O	SER A 471	-41.761	30.899	16.067	1.00	25.43	O
ATOM	7177	N	VAL A 472	-42.653	29.381	17.482	1.00	24.95	N
ATOM	7178	CA	VAL A 472	-43.097	28.451	16.438	1.00	24.65	C
ATOM	7180	CB	VAL A 472	-43.485	27.055	17.024	1.00	24.57	C
ATOM	7182	CG1	VAL A 472	-44.257	26.222	16.016	1.00	23.23	C
ATOM	7186	CG2	VAL A 472	-42.232	26.309	17.494	1.00	23.83	C
ATOM	7190	C	VAL A 472	-44.252	29.039	15.631	1.00	24.71	C
ATOM	7191	O	VAL A 472	-44.369	28.776	14.440	1.00	24.63	O
ATOM	7193	N	MET A 473	-45.073	29.866	16.269	1.00	24.79	N
ATOM	7194	CA	MET A 473	-46.143	30.591	15.565	1.00	24.94	C
ATOM	7196	CB	MET A 473	-47.059	31.266	16.576	1.00	25.11	C
ATOM	7199	CG	MET A 473	-47.683	30.335	17.560	1.00	25.11	C
ATOM	7202	SD	MET A 473	-48.967	29.391	16.780	1.00	25.35	S
ATOM	7203	CE	MET A 473	-50.064	29.163	18.189	1.00	26.11	C

ATOM	7207	C	MET A 473	-45.613	31.672	14.605	1.00	24.86	C
ATOM	7208	O	MET A 473	-46.132	31.849	13.513	1.00	24.57	O
ATOM	7210	N	ASN A 474	-44.589	32.405	15.032	1.00	24.89	N
ATOM	7211	CA	ASN A 474	-43.951	33.383	14.166	1.00	25.00	C
ATOM	7213	CB	ASN A 474	-43.009	34.278	14.966	1.00	25.12	C
ATOM	7216	CG	ASN A 474	-43.700	34.966	16.127	1.00	25.77	C
ATOM	7217	OD1	ASN A 474	-43.058	35.312	17.114	1.00	26.96	O
ATOM	7218	ND2	ASN A 474	-45.015	35.163	16.021	1.00	26.31	N
ATOM	7221	C	ASN A 474	-43.195	32.707	13.024	1.00	24.84	C
ATOM	7222	O	ASN A 474	-42.979	33.313	11.973	1.00	25.07	O
ATOM	7224	N	LEU A 475	-42.792	31.455	13.227	1.00	24.46	N
ATOM	7225	CA	LEU A 475	-42.218	30.671	12.142	1.00	24.06	C
ATOM	7227	CB	LEU A 475	-41.596	29.367	12.653	1.00	24.27	C
ATOM	7230	CG	LEU A 475	-40.307	29.072	11.886	1.00	25.08	C
ATOM	7232	CD1	LEU A 475	-39.207	30.000	12.419	1.00	25.89	C
ATOM	7236	CD2	LEU A 475	-39.887	27.622	11.985	1.00	25.33	C
ATOM	7240	C	LEU A 475	-43.278	30.365	11.085	1.00	23.26	C
ATOM	7241	O	LEU A 475	-43.018	30.479	9.884	1.00	22.90	O
ATOM	7243	N	ILE A 476	-44.471	29.978	11.534	1.00	22.53	N
ATOM	7244	CA	ILE A 476	-45.544	29.632	10.606	1.00	21.93	C
ATOM	7246	CB	ILE A 476	-46.773	29.017	11.308	1.00	21.58	C
ATOM	7248	CG1	ILE A 476	-46.499	27.550	11.664	1.00	21.32	C
ATOM	7251	CD1	ILE A 476	-47.552	26.902	12.598	1.00	20.33	C
ATOM	7255	CG2	ILE A 476	-47.994	29.104	10.428	1.00	20.46	C
ATOM	7259	C	ILE A 476	-45.926	30.887	9.853	1.00	22.07	C
ATOM	7260	O	ILE A 476	-46.007	30.869	8.626	1.00	22.47	O
ATOM	7262	N	ASP A 477	-46.122	31.982	10.584	1.00	21.81	N
ATOM	7263	CA	ASP A 477	-46.483	33.245	9.970	1.00	21.56	C
ATOM	7265	CB	ASP A 477	-46.643	34.337	11.032	1.00	21.89	C
ATOM	7268	CG	ASP A 477	-47.962	34.217	11.817	1.00	23.39	C
ATOM	7269	OD1	ASP A 477	-48.854	33.441	11.375	1.00	24.41	O
ATOM	7270	OD2	ASP A 477	-48.103	34.902	12.875	1.00	23.80	O
ATOM	7271	C	ASP A 477	-45.425	33.626	8.947	1.00	21.06	C
ATOM	7272	O	ASP A 477	-45.759	33.913	7.795	1.00	20.78	O
ATOM	7274	N	GLU A 478	-44.156	33.591	9.361	1.00	20.63	N
ATOM	7275	CA	GLU A 478	-43.027	33.919	8.467	1.00	20.45	C
ATOM	7277	CB	GLU A 478	-41.680	33.793	9.200	1.00	20.66	C
ATOM	7280	CG	GLU A 478	-40.422	34.044	8.339	1.00	22.27	C
ATOM	7283	CD	GLU A 478	-39.107	33.579	9.018	1.00	24.72	C
ATOM	7284	OE1	GLU A 478	-38.553	34.367	9.817	1.00	26.19	O
ATOM	7285	OE2	GLU A 478	-38.619	32.442	8.745	1.00	25.66	O
ATOM	7286	C	GLU A 478	-43.041	33.042	7.217	1.00	19.54	C
ATOM	7287	O	GLU A 478	-42.879	33.544	6.102	1.00	19.02	O
ATOM	7289	N	THR A 479	-43.250	31.743	7.413	1.00	18.76	N
ATOM	7290	CA	THR A 479	-43.358	30.808	6.294	1.00	18.45	C
ATOM	7292	CB	THR A 479	-43.503	29.339	6.774	1.00	18.26	C
ATOM	7294	OG1	THR A 479	-42.334	28.955	7.494	1.00	17.30	O
ATOM	7296	CG2	THR A 479	-43.669	28.394	5.592	1.00	17.93	C
ATOM	7300	C	THR A 479	-44.504	31.173	5.314	1.00	18.39	C
ATOM	7301	O	THR A 479	-44.280	31.200	4.103	1.00	18.58	O
ATOM	7303	N	TRP A 480	-45.709	31.457	5.818	1.00	17.85	N
ATOM	7304	CA	TRP A 480	-46.801	31.909	4.942	1.00	17.44	C
ATOM	7306	CB	TRP A 480	-48.086	32.225	5.728	1.00	17.67	C
ATOM	7309	CG	TRP A 480	-48.969	31.042	5.888	1.00	17.09	C
ATOM	7310	CD1	TRP A 480	-48.996	30.195	6.936	1.00	16.76	C
ATOM	7312	NE1	TRP A 480	-49.914	29.216	6.724	1.00	17.08	N
ATOM	7314	CE2	TRP A 480	-50.501	29.410	5.507	1.00	17.86	C
ATOM	7315	CD2	TRP A 480	-49.930	30.559	4.954	1.00	17.15	C
ATOM	7316	CE3	TRP A 480	-50.356	30.980	3.696	1.00	18.31	C
ATOM	7318	CZ3	TRP A 480	-51.337	30.244	3.039	1.00	19.10	C

ATOM	7320	CH2 TRP A 480	-51.897	29.107	3.622	1.00	18.83	C
ATOM	7322	CZ2 TRP A 480	-51.493	28.673	4.855	1.00	18.86	C
ATOM	7324	C TRP A 480	-46.426	33.124	4.097	1.00	17.10	C
ATOM	7325	O TRP A 480	-46.824	33.215	2.943	1.00	17.22	O
ATOM	7327	N LYS A 481	-45.675	34.062	4.659	1.00	16.60	N
ATOM	7328	CA LYS A 481	-45.275	35.215	3.886	1.00	16.23	C
ATOM	7330	CB LYS A 481	-44.566	36.243	4.747	1.00	16.42	C
ATOM	7333	CG LYS A 481	-45.417	36.954	5.769	1.00	16.42	C
ATOM	7336	CD LYS A 481	-44.555	37.979	6.505	1.00	16.40	C
ATOM	7339	CE LYS A 481	-45.199	38.504	7.782	1.00	16.66	C
ATOM	7342	NZ LYS A 481	-44.184	38.717	8.846	1.00	16.72	N
ATOM	7346	C LYS A 481	-44.357	34.782	2.767	1.00	15.94	C
ATOM	7347	O LYS A 481	-44.451	35.291	1.674	1.00	16.09	O
ATOM	7349	N LYS A 482	-43.461	33.849	3.034	1.00	15.98	N
ATOM	7350	CA LYS A 482	-42.559	33.369	1.996	1.00	16.30	C
ATOM	7352	CB LYS A 482	-41.409	32.563	2.594	1.00	16.51	C
ATOM	7355	CG LYS A 482	-40.354	33.445	3.286	1.00	17.46	C
ATOM	7358	CD LYS A 482	-39.486	32.647	4.241	1.00	18.48	C
ATOM	7361	CE LYS A 482	-38.625	33.548	5.096	1.00	18.94	C
ATOM	7364	NZ LYS A 482	-37.865	32.773	6.129	1.00	20.40	N
ATOM	7368	C LYS A 482	-43.296	32.561	.938	1.00	16.40	C
ATOM	7369	O LYS A 482	-42.941	32.614	-.229	1.00	16.85	O
ATOM	7371	N MET A 483	-44.328	31.824	1.332	1.00	16.57	N
ATOM	7372	CA MET A 483	-45.149	31.100	.366	1.00	16.67	C
ATOM	7374	CB MET A 483	-46.128	30.156	1.057	1.00	16.55	C
ATOM	7377	CG MET A 483	-45.496	28.923	1.675	1.00	16.50	C
ATOM	7380	SD MET A 483	-46.715	27.684	2.194	1.00	17.39	S
ATOM	7381	CE MET A 483	-47.937	28.704	3.025	1.00	17.27	C
ATOM	7385	C MET A 483	-45.928	32.069	-.495	1.00	17.06	C
ATOM	7386	O MET A 483	-46.027	31.867	-1.688	1.00	17.17	O
ATOM	7388	N ASN A 484	-46.483	33.110	.118	1.00	17.72	N
ATOM	7389	CA ASN A 484	-47.305	34.093	-.587	1.00	18.35	C
ATOM	7391	CB ASN A 484	-47.861	35.129	.397	1.00	18.26	C
ATOM	7394	CG ASN A 484	-48.994	34.601	1.248	1.00	17.35	C
ATOM	7395	OD1 ASN A 484	-49.597	33.577	.934	1.00	17.54	O
ATOM	7396	ND2 ASN A 484	-49.304	35.317	2.331	1.00	14.10	N
ATOM	7399	C ASN A 484	-46.541	34.839	-1.677	1.00	19.51	C
ATOM	7400	O ASN A 484	-47.123	35.250	-2.685	1.00	19.49	O
ATOM	7402	N LYS A 485	-45.243	35.035	-1.463	1.00	20.86	N
ATOM	7403	CA LYS A 485	-44.394	35.681	-2.456	1.00	22.12	C
ATOM	7405	CB LYS A 485	-43.070	36.092	-1.823	1.00	22.08	C
ATOM	7408	CG LYS A 485	-42.229	37.060	-2.641	1.00	23.09	C
ATOM	7411	CD LYS A 485	-40.783	37.072	-2.109	1.00	25.38	C
ATOM	7414	CE LYS A 485	-40.098	38.460	-2.135	1.00	26.12	C
ATOM	7417	NZ LYS A 485	-38.943	38.555	-3.086	1.00	26.61	N
ATOM	7421	C LYS A 485	-44.164	34.743	-3.656	1.00	23.36	C
ATOM	7422	O LYS A 485	-44.120	35.191	-4.802	1.00	23.00	O
ATOM	7424	N GLU A 486	-44.027	33.442	-3.400	1.00	24.99	N
ATOM	7425	CA GLU A 486	-43.902	32.479	-4.491	1.00	26.35	C
ATOM	7427	CB GLU A 486	-43.627	31.061	-3.978	1.00	26.55	C
ATOM	7430	CG GLU A 486	-42.985	30.130	-5.033	1.00	28.29	C
ATOM	7433	CD GLU A 486	-41.466	30.324	-5.191	1.00	31.11	C
ATOM	7434	OE1 GLU A 486	-40.981	30.466	-6.341	1.00	32.07	O
ATOM	7435	OE2 GLU A 486	-40.744	30.335	-4.159	1.00	33.04	O
ATOM	7436	C GLU A 486	-45.157	32.501	-5.362	1.00	27.13	C
ATOM	7437	O GLU A 486	-45.064	32.618	-6.583	1.00	27.57	O
ATOM	7439	N LYS A 487	-46.325	32.422	-4.739	1.00	28.12	N
ATOM	7440	CA LYS A 487	-47.582	32.378	-5.485	1.00	29.12	C
ATOM	7442	CB LYS A 487	-48.788	32.260	-4.542	1.00	29.21	C
ATOM	7445	CG LYS A 487	-50.167	32.296	-5.235	1.00	28.57	C

ATOM	7448	CD	LYS	A 487	-50.403	31.040	-6.036	1.00	27.91	C
ATOM	7451	CE	LYS	A 487	-51.615	31.141	-6.936	1.00	28.38	C
ATOM	7454	NZ	LYS	A 487	-52.874	31.460	-6.211	1.00	28.03	N
ATOM	7458	C	LYS	A 487	-47.756	33.603	-6.352	1.00	30.18	C
ATOM	7459	O	LYS	A 487	-48.308	33.514	-7.446	1.00	30.21	O
ATOM	7461	N	LEU	A 488	-47.295	34.743	-5.844	1.00	31.64	N
ATOM	7462	CA	LEU	A 488	-47.422	36.026	-6.532	1.00	32.65	C
ATOM	7464	CB	LEU	A 488	-47.426	37.144	-5.494	1.00	32.43	C
ATOM	7467	CG	LEU	A 488	-48.091	38.446	-5.907	1.00	32.50	C
ATOM	7469	CD1	LEU	A 488	-49.593	38.288	-5.943	1.00	32.66	C
ATOM	7473	CD2	LEU	A 488	-47.707	39.554	-4.943	1.00	33.38	C
ATOM	7477	C	LEU	A 488	-46.286	36.242	-7.542	1.00	33.99	C
ATOM	7478	O	LEU	A 488	-46.536	36.539	-8.705	1.00	33.72	O
ATOM	7480	N	GLY	A 489	-45.046	36.045	-7.091	1.00	35.85	N
ATOM	7481	CA	GLY	A 489	-43.845	36.361	-7.874	1.00	37.37	C
ATOM	7484	C	GLY	A 489	-43.278	35.221	-8.709	1.00	38.76	C
ATOM	7485	O	GLY	A 489	-42.414	34.467	-8.250	1.00	39.06	O
ATOM	7487	N	GLY	A 490	-43.760	35.119	-9.947	1.00	40.26	N
ATOM	7488	CA	GLY	A 490	-43.285	34.130	-10.927	1.00	41.03	C
ATOM	7491	C	GLY	A 490	-42.193	33.194	-10.436	1.00	41.61	C
ATOM	7492	O	GLY	A 490	-41.010	33.547	-10.423	1.00	41.62	O
ATOM	7494	N	SER	A 491	-42.590	31.996	-10.023	1.00	42.17	N
ATOM	7495	CA	SER	A 491	-41.619	30.962	-9.662	1.00	42.54	C
ATOM	7497	CB	SER	A 491	-42.279	29.909	-8.741	1.00	42.63	C
ATOM	7500	OG	SER	A 491	-43.451	29.345	-9.294	1.00	42.44	O
ATOM	7502	C	SER	A 491	-41.032	30.349	-10.952	1.00	42.50	C
ATOM	7503	O	SER	A 491	-41.117	30.961	-12.030	1.00	42.57	O
ATOM	7505	N	LEU	A 492	-40.405	29.177	-10.843	1.00	42.18	N
ATOM	7506	CA	LEU	A 492	-40.185	28.333	-12.021	1.00	41.96	C
ATOM	7508	CB	LEU	A 492	-39.110	27.270	-11.782	1.00	42.39	C
ATOM	7511	CG	LEU	A 492	-37.695	27.671	-11.353	1.00	44.21	C
ATOM	7513	CD1	LEU	A 492	-36.783	26.438	-11.498	1.00	45.26	C
ATOM	7517	CD2	LEU	A 492	-37.120	28.896	-12.129	1.00	45.51	C
ATOM	7521	C	LEU	A 492	-41.479	27.610	-12.356	1.00	40.88	C
ATOM	7522	O	LEU	A 492	-41.713	27.244	-13.504	1.00	40.78	O
ATOM	7524	N	PHE	A 493	-42.308	27.409	-11.339	1.00	39.61	N
ATOM	7525	CA	PHE	A 493	-43.500	26.595	-11.460	1.00	38.85	C
ATOM	7527	CB	PHE	A 493	-43.783	25.911	-10.124	1.00	38.62	C
ATOM	7530	CG	PHE	A 493	-42.725	24.937	-9.713	1.00	37.16	C
ATOM	7531	CD1	PHE	A 493	-42.901	23.587	-9.916	1.00	35.26	C
ATOM	7533	CE1	PHE	A 493	-41.934	22.698	-9.546	1.00	34.86	C
ATOM	7535	CZ	PHE	A 493	-40.769	23.143	-8.966	1.00	35.09	C
ATOM	7537	CE2	PHE	A 493	-40.574	24.482	-8.759	1.00	35.45	C
ATOM	7539	CD2	PHE	A 493	-41.549	25.374	-9.130	1.00	36.36	C
ATOM	7541	C	PHE	A 493	-44.714	27.412	-11.890	1.00	38.62	C
ATOM	7542	O	PHE	A 493	-44.756	28.621	-11.702	1.00	38.86	O
ATOM	7544	N	ALA	A 494	-45.698	26.731	-12.469	1.00	38.32	N
ATOM	7545	CA	ALA	A 494	-46.977	27.337	-12.843	1.00	38.09	C
ATOM	7547	CB	ALA	A 494	-47.658	26.486	-13.906	1.00	38.17	C
ATOM	7551	C	ALA	A 494	-47.891	27.479	-11.623	1.00	37.76	C
ATOM	7552	O	ALA	A 494	-48.039	26.541	-10.845	1.00	38.15	O
ATOM	7554	N	LYS	A 495	-48.541	28.629	-11.484	1.00	37.15	N
ATOM	7555	CA	LYS	A 495	-49.321	28.947	-10.270	1.00	36.53	C
ATOM	7557	CB	LYS	A 495	-50.116	30.258	-10.477	1.00	36.90	C
ATOM	7560	CG	LYS	A 495	-49.235	31.524	-10.582	1.00	37.72	C
ATOM	7563	CD	LYS	A 495	-50.061	32.826	-10.585	1.00	38.87	C
ATOM	7566	CE	LYS	A 495	-49.140	34.065	-10.665	1.00	39.78	C
ATOM	7569	NZ	LYS	A 495	-49.786	35.369	-10.282	1.00	39.78	N
ATOM	7573	C	LYS	A 495	-50.238	27.809	-9.724	1.00	35.25	C
ATOM	7574	O	LYS	A 495	-50.261	27.563	-8.523	1.00	34.91	O

ATOM	7576	N	PRO A 496	-50.978	27.107	-10.599	1.00	33.76	N
ATOM	7577	CA	PRO A 496	-51.846	26.015	-10.162	1.00	32.74	C
ATOM	7579	CB	PRO A 496	-52.222	25.340	-11.474	1.00	32.99	C
ATOM	7582	CG	PRO A 496	-52.337	26.470	-12.385	1.00	33.75	C
ATOM	7585	CD	PRO A 496	-51.196	27.400	-12.022	1.00	33.83	C
ATOM	7588	C	PRO A 496	-51.194	24.998	-9.267	1.00	31.24	C
ATOM	7589	O	PRO A 496	-51.822	24.522	-8.331	1.00	31.54	O
ATOM	7590	N	PHE A 497	-49.954	24.642	-9.566	1.00	29.40	N
ATOM	7591	CA	PHE A 497	-49.229	23.726	-8.708	1.00	27.76	C
ATOM	7593	CB	PHE A 497	-48.162	22.948	-9.474	1.00	27.62	C
ATOM	7596	CG	PHE A 497	-47.351	22.040	-8.597	1.00	26.33	C
ATOM	7597	CD1	PHE A 497	-47.944	20.969	-7.974	1.00	24.96	C
ATOM	7599	CE1	PHE A 497	-47.218	20.151	-7.156	1.00	25.04	C
ATOM	7601	CZ	PHE A 497	-45.875	20.396	-6.940	1.00	24.79	C
ATOM	7603	CE2	PHE A 497	-45.274	21.457	-7.549	1.00	24.73	C
ATOM	7605	CD2	PHE A 497	-46.010	22.282	-8.366	1.00	25.48	C
ATOM	7607	C	PHE A 497	-48.592	24.407	-7.494	1.00	26.59	C
ATOM	7608	O	PHE A 497	-48.361	23.738	-6.500	1.00	26.77	O
ATOM	7610	N	VAL A 498	-48.291	25.704	-7.535	1.00	24.83	N
ATOM	7611	CA	VAL A 498	-47.869	26.338	-6.289	1.00	23.86	C
ATOM	7613	CB	VAL A 498	-47.181	27.722	-6.455	1.00	23.47	C
ATOM	7615	CG1	VAL A 498	-48.163	28.755	-6.759	1.00	24.57	C
ATOM	7619	CG2	VAL A 498	-46.151	27.690	-7.551	1.00	23.54	C
ATOM	7623	C	VAL A 498	-49.094	26.396	-5.355	1.00	22.97	C
ATOM	7624	O	VAL A 498	-48.978	26.130	-4.161	1.00	23.39	O
ATOM	7626	N	GLU A 499	-50.268	26.687	-5.903	1.00	21.68	N
ATOM	7627	CA	GLU A 499	-51.473	26.719	-5.106	1.00	20.84	C
ATOM	7629	CB	GLU A 499	-52.677	27.136	-5.930	1.00	20.99	C
ATOM	7632	CG	GLU A 499	-53.957	27.352	-5.099	1.00	21.25	C
ATOM	7635	CD	GLU A 499	-53.930	28.630	-4.265	1.00	21.77	C
ATOM	7636	OE1	GLU A 499	-52.982	29.448	-4.383	1.00	20.74	O
ATOM	7637	OE2	GLU A 499	-54.876	28.812	-3.476	1.00	22.76	O
ATOM	7638	C	GLU A 499	-51.766	25.371	-4.497	1.00	20.35	C
ATOM	7639	O	GLU A 499	-52.137	25.308	-3.318	1.00	20.78	O
ATOM	7641	N	THR A 500	-51.629	24.286	-5.264	1.00	19.40	N
ATOM	7642	CA	THR A 500	-51.894	22.970	-4.667	1.00	18.94	C
ATOM	7644	CB	THR A 500	-51.993	21.782	-5.683	1.00	18.96	C
ATOM	7646	OG1	THR A 500	-50.708	21.475	-6.202	1.00	19.44	O
ATOM	7648	CG2	THR A 500	-52.991	22.071	-6.835	1.00	18.75	C
ATOM	7652	C	THR A 500	-50.895	22.674	-3.526	1.00	18.20	C
ATOM	7653	O	THR A 500	-51.280	22.100	-2.519	1.00	17.68	O
ATOM	7655	N	ALA A 501	-49.646	23.119	-3.666	1.00	17.49	N
ATOM	7656	CA	ALA A 501	-48.663	23.021	-2.593	1.00	17.18	C
ATOM	7658	CB	ALA A 501	-47.348	23.563	-3.051	1.00	16.98	C
ATOM	7662	C	ALA A 501	-49.120	23.770	-1.338	1.00	17.38	C
ATOM	7663	O	ALA A 501	-49.098	23.219	-.230	1.00	17.45	O
ATOM	7665	N	ILE A 502	-49.517	25.031	-1.508	1.00	17.29	N
ATOM	7666	CA	ILE A 502	-49.971	25.857	-.386	1.00	17.07	C
ATOM	7668	CB	ILE A 502	-50.353	27.304	-.846	1.00	17.09	C
ATOM	7670	CG1	ILE A 502	-49.092	28.069	-1.286	1.00	17.17	C
ATOM	7673	CD1	ILE A 502	-49.345	29.345	-2.116	1.00	16.25	C
ATOM	7677	CG2	ILE A 502	-51.110	28.076	.265	1.00	16.14	C
ATOM	7681	C	ILE A 502	-51.155	25.168	.285	1.00	17.24	C
ATOM	7682	O	ILE A 502	-51.265	25.163	1.516	1.00	17.03	O
ATOM	7684	N	ASN A 503	-52.022	24.559	-.522	1.00	17.32	N
ATOM	7685	CA	ASN A 503	-53.167	23.823	.022	1.00	17.75	C
ATOM	7687	CB	ASN A 503	-53.986	23.221	-1.121	1.00	17.70	C
ATOM	7690	CG	ASN A 503	-54.760	24.261	-1.888	1.00	18.41	C
ATOM	7691	OD1	ASN A 503	-55.058	25.334	-1.382	1.00	18.94	O
ATOM	7692	ND2	ASN A 503	-55.107	23.939	-3.119	1.00	20.42	N

ATOM	7695	C	ASN A 503	-52.803	22.727	1.075	1.00	17.74	C
ATOM	7696	O	ASN A 503	-53.619	22.387	1.949	1.00	17.96	O
ATOM	7698	N	LEU A 504	-51.589	22.185	.993	1.00	17.23	N
ATOM	7699	CA	LEU A 504	-51.105	21.249	1.994	1.00	17.10	C
ATOM	7701	CB	LEU A 504	-49.745	20.686	1.583	1.00	17.28	C
ATOM	7704	CG	LEU A 504	-49.213	19.528	2.426	1.00	17.29	C
ATOM	7706	CD1	LEU A 504	-49.570	18.196	1.762	1.00	17.22	C
ATOM	7710	CD2	LEU A 504	-47.708	19.686	2.616	1.00	16.83	C
ATOM	7714	C	LEU A 504	-50.971	21.939	3.347	1.00	16.97	C
ATOM	7715	O	LEU A 504	-51.237	21.330	4.379	1.00	16.84	O
ATOM	7717	N	ALA A 505	-50.535	23.201	3.330	1.00	16.90	N
ATOM	7718	CA	ALA A 505	-50.433	24.015	4.540	1.00	16.68	C
ATOM	7720	CB	ALA A 505	-49.739	25.309	4.243	1.00	16.54	C
ATOM	7724	C	ALA A 505	-51.826	24.281	5.074	1.00	16.77	C
ATOM	7725	O	ALA A 505	-52.087	24.123	6.266	1.00	16.80	O
ATOM	7727	N	ARG A 506	-52.726	24.664	4.172	1.00	16.82	N
ATOM	7728	CA	ARG A 506	-54.128	24.875	4.521	1.00	16.84	C
ATOM	7730	CB	ARG A 506	-54.944	25.274	3.286	1.00	16.86	C
ATOM	7733	CG	ARG A 506	-54.649	26.661	2.795	1.00	16.27	C
ATOM	7736	CD	ARG A 506	-55.586	27.090	1.726	1.00	15.26	C
ATOM	7739	NE	ARG A 506	-55.240	28.436	1.273	1.00	15.31	N
ATOM	7741	CZ	ARG A 506	-54.744	28.758	.082	1.00	14.28	C
ATOM	7742	NH1	ARG A 506	-54.519	27.856	-.856	1.00	14.12	N
ATOM	7745	NH2	ARG A 506	-54.471	30.018	-.176	1.00	15.94	N
ATOM	7748	C	ARG A 506	-54.732	23.631	5.122	1.00	16.96	C
ATOM	7749	O	ARG A 506	-55.480	23.707	6.089	1.00	16.71	O
ATOM	7751	N	GLN A 507	-54.415	22.481	4.542	1.00	17.37	N
ATOM	7752	CA	GLN A 507	-54.962	21.233	5.048	1.00	17.87	C
ATOM	7754	CB	GLN A 507	-54.712	20.075	4.087	1.00	17.80	C
ATOM	7757	CG	GLN A 507	-55.293	18.740	4.571	1.00	17.24	C
ATOM	7760	CD	GLN A 507	-56.777	18.805	4.860	1.00	16.23	C
ATOM	7761	OE1	GLN A 507	-57.506	19.585	4.250	1.00	16.66	O
ATOM	7762	NE2	GLN A 507	-57.233	17.980	5.790	1.00	14.96	N
ATOM	7765	C	GLN A 507	-54.401	20.890	6.423	1.00	18.28	C
ATOM	7766	O	GLN A 507	-55.148	20.419	7.287	1.00	18.53	O
ATOM	7768	N	SER A 508	-53.096	21.113	6.609	1.00	18.50	N
ATOM	7769	CA	SER A 508	-52.435	20.924	7.908	1.00	18.54	C
ATOM	7771	CB	SER A 508	-50.979	21.385	7.856	1.00	18.40	C
ATOM	7774	OG	SER A 508	-50.259	20.694	6.857	1.00	18.51	O
ATOM	7776	C	SER A 508	-53.147	21.752	8.942	1.00	18.78	C
ATOM	7777	O	SER A 508	-53.535	21.278	10.006	1.00	18.33	O
ATOM	7779	N	HIS A 509	-53.324	23.014	8.599	1.00	19.37	N
ATOM	7780	CA	HIS A 509	-54.014	23.922	9.466	1.00	19.92	C
ATOM	7782	CB	HIS A 509	-54.126	25.295	8.818	1.00	20.07	C
ATOM	7785	CG	HIS A 509	-53.000	26.188	9.174	1.00	19.82	C
ATOM	7786	ND1	HIS A 509	-52.079	26.631	8.255	1.00	19.73	N
ATOM	7788	CE1	HIS A 509	-51.187	27.383	8.866	1.00	20.25	C
ATOM	7790	NE2	HIS A 509	-51.485	27.426	10.149	1.00	21.66	N
ATOM	7792	CD2	HIS A 509	-52.615	26.682	10.368	1.00	20.98	C
ATOM	7794	C	HIS A 509	-55.389	23.448	9.839	1.00	20.41	C
ATOM	7795	O	HIS A 509	-55.821	23.694	10.934	1.00	20.22	O
ATOM	7797	N	CYS A 510	-56.093	22.808	8.918	1.00	21.33	N
ATOM	7798	CA	CYS A 510	-57.470	22.396	9.189	1.00	22.03	C
ATOM	7800	CB	CYS A 510	-58.278	22.383	7.894	1.00	21.90	C
ATOM	7803	SG	CYS A 510	-58.522	24.055	7.270	1.00	21.93	S
ATOM	7805	C	CYS A 510	-57.527	21.054	9.908	1.00	22.60	C
ATOM	7806	O	CYS A 510	-58.423	20.822	10.706	1.00	22.23	O
ATOM	7808	N	THR A 511	-56.553	20.194	9.627	1.00	23.67	N
ATOM	7809	CA	THR A 511	-56.421	18.920	10.305	1.00	24.81	C
ATOM	7811	CB	THR A 511	-55.437	18.000	9.565	1.00	24.70	C

ATOM	7813	OG1 THR A 511	-56.071	17.485	8.394	1.00	24.45	O
ATOM	7815	CG2 THR A 511	-55.003	16.844	10.448	1.00	24.23	C
ATOM	7819	C THR A 511	-55.965	19.086	11.756	1.00	26.22	C
ATOM	7820	O THR A 511	-56.685	18.717	12.679	1.00	26.12	O
ATOM	7822	N TYR A 512	-54.775	19.643	11.959	1.00	28.08	N
ATOM	7823	CA TYR A 512	-54.177	19.657	13.298	1.00	29.63	C
ATOM	7825	CB TYR A 512	-52.663	19.783	13.238	1.00	29.67	C
ATOM	7828	CG TYR A 512	-52.117	18.579	12.560	1.00	29.56	C
ATOM	7829	CD1 TYR A 512	-51.965	17.395	13.246	1.00	29.66	C
ATOM	7831	CE1 TYR A 512	-51.507	16.273	12.614	1.00	30.94	C
ATOM	7833	CZ TYR A 512	-51.217	16.326	11.263	1.00	32.04	C
ATOM	7834	OH TYR A 512	-50.757	15.209	10.608	1.00	33.78	O
ATOM	7836	CE2 TYR A 512	-51.390	17.492	10.557	1.00	31.29	C
ATOM	7838	CD2 TYR A 512	-51.847	18.600	11.205	1.00	30.53	C
ATOM	7840	C TYR A 512	-54.812	20.676	14.206	1.00	31.11	C
ATOM	7841	O TYR A 512	-55.494	20.277	15.139	1.00	31.34	O
ATOM	7843	N HIS A 513	-54.583	21.969	13.966	1.00	32.89	N
ATOM	7844	CA HIS A 513	-55.505	23.023	14.442	1.00	34.48	C
ATOM	7846	CB HIS A 513	-56.372	23.461	13.225	1.00	35.13	C
ATOM	7849	CG HIS A 513	-57.704	24.108	13.529	1.00	37.10	C
ATOM	7850	ND1 HIS A 513	-58.265	25.040	12.678	1.00	38.69	N
ATOM	7852	CE1 HIS A 513	-59.442	25.416	13.151	1.00	39.23	C
ATOM	7854	NE2 HIS A 513	-59.678	24.751	14.268	1.00	38.18	N
ATOM	7856	CD2 HIS A 513	-58.615	23.915	14.519	1.00	38.17	C
ATOM	7858	C HIS A 513	-56.337	22.503	15.610	1.00	35.00	C
ATOM	7859	O HIS A 513	-56.302	23.074	16.712	1.00	35.48	O
ATOM	7861	N ASN A 514	-57.079	21.419	15.351	1.00	35.23	N
ATOM	7862	CA ASN A 514	-57.934	20.768	16.347	1.00	35.46	C
ATOM	7864	CB ASN A 514	-58.342	19.343	15.903	1.00	35.39	C
ATOM	7867	CG ASN A 514	-59.383	19.345	14.786	1.00	34.01	C
ATOM	7868	OD1 ASN A 514	-59.336	20.185	13.896	1.00	33.36	O
ATOM	7869	ND2 ASN A 514	-60.315	18.406	14.833	1.00	31.13	N
ATOM	7872	C ASN A 514	-57.464	20.751	17.818	1.00	36.15	C
ATOM	7873	O ASN A 514	-56.397	20.211	18.153	1.00	36.13	O
ATOM	7875	N GLY A 515	-58.277	21.430	18.640	1.00	36.95	N
ATOM	7876	CA GLY A 515	-58.427	21.235	20.079	1.00	37.43	C
ATOM	7879	C GLY A 515	-59.724	20.470	20.381	1.00	38.18	C
ATOM	7880	O GLY A 515	-59.636	19.418	21.021	1.00	38.57	O
ATOM	7882	N ASP A 516	-60.930	20.921	19.969	1.00	38.73	N
ATOM	7883	CA ASP A 516	-61.249	22.194	19.293	1.00	39.26	C
ATOM	7885	CB ASP A 516	-61.607	21.904	17.845	1.00	39.51	C
ATOM	7888	CG ASP A 516	-60.523	22.291	16.894	1.00	42.32	C
ATOM	7889	OD1 ASP A 516	-59.613	23.084	17.286	1.00	45.56	O
ATOM	7890	OD2 ASP A 516	-60.578	21.797	15.740	1.00	45.52	O
ATOM	7891	C ASP A 516	-62.466	22.913	19.884	1.00	39.38	C
ATOM	7892	O ASP A 516	-63.032	22.454	20.880	1.00	39.74	O
ATOM	7894	N ALA A 517	-62.850	24.038	19.259	1.00	39.40	N
ATOM	7895	CA ALA A 517	-64.140	24.752	19.481	1.00	39.50	C
ATOM	7897	CB ALA A 517	-65.117	24.403	18.333	1.00	39.27	C
ATOM	7901	C ALA A 517	-64.837	24.566	20.863	1.00	39.90	C
ATOM	7902	O ALA A 517	-64.173	24.428	21.899	1.00	39.89	O
ATOM	7904	N HIS A 518	-66.175	24.617	20.876	1.00	40.33	N
ATOM	7905	CA HIS A 518	-66.973	24.112	22.016	1.00	40.73	C
ATOM	7907	CB HIS A 518	-68.091	25.092	22.396	1.00	41.22	C
ATOM	7910	CG HIS A 518	-67.578	26.419	22.866	1.00	43.58	C
ATOM	7911	ND1 HIS A 518	-67.426	26.725	24.207	1.00	45.88	N
ATOM	7913	CE1 HIS A 518	-66.938	27.950	24.320	1.00	46.42	C
ATOM	7915	NE2 HIS A 518	-66.760	28.446	23.102	1.00	46.25	N
ATOM	7917	CD2 HIS A 518	-67.145	27.506	22.174	1.00	45.24	C
ATOM	7919	C HIS A 518	-67.526	22.725	21.655	1.00	40.17	C

ATOM	7920	O	HIS A 518	-68.738	22.501	21.571	1.00	40.02	O
ATOM	7922	N	THR A 519	-66.581	21.812	21.448	1.00	39.64	N
ATOM	7923	CA	THR A 519	-66.808	20.466	20.918	1.00	39.16	C
ATOM	7925	CB	THR A 519	-67.414	20.469	19.453	1.00	39.16	C
ATOM	7927	OG1	THR A 519	-66.847	21.529	18.665	1.00	39.04	O
ATOM	7929	CG2	THR A 519	-68.942	20.633	19.480	1.00	38.82	C
ATOM	7933	C	THR A 519	-65.428	19.771	20.978	1.00	38.85	C
ATOM	7934	O	THR A 519	-64.391	20.433	20.878	1.00	38.66	O
ATOM	7936	N	SER A 520	-65.407	18.455	21.170	1.00	38.37	N
ATOM	7937	CA	SER A 520	-64.160	17.745	21.500	1.00	38.08	C
ATOM	7939	CB	SER A 520	-64.514	16.320	21.931	1.00	38.06	C
ATOM	7942	OG	SER A 520	-64.560	15.475	20.805	1.00	38.66	O
ATOM	7944	C	SER A 520	-63.136	17.771	20.327	1.00	37.86	C
ATOM	7945	O	SER A 520	-63.420	18.369	19.289	1.00	37.50	O
ATOM	7947	N	PRO A 521	-61.948	17.122	20.482	1.00	37.98	N
ATOM	7948	CA	PRO A 521	-60.944	17.170	19.379	1.00	37.92	C
ATOM	7950	CB	PRO A 521	-59.674	16.527	19.987	1.00	37.87	C
ATOM	7953	CG	PRO A 521	-60.091	15.919	21.330	1.00	38.28	C
ATOM	7956	CD	PRO A 521	-61.577	16.137	21.526	1.00	38.13	C
ATOM	7959	C	PRO A 521	-61.432	16.446	18.110	1.00	37.96	C
ATOM	7960	O	PRO A 521	-61.769	17.118	17.140	1.00	38.09	O
ATOM	7961	N	ASP A 522	-61.483	15.108	18.107	1.00	38.07	N
ATOM	7962	CA	ASP A 522	-62.351	14.393	17.161	1.00	38.14	C
ATOM	7964	CB	ASP A 522	-62.132	12.877	17.204	1.00	38.31	C
ATOM	7967	CG	ASP A 522	-60.769	12.456	16.640	1.00	39.45	C
ATOM	7968	OD1	ASP A 522	-60.160	13.200	15.827	1.00	39.72	O
ATOM	7969	OD2	ASP A 522	-60.302	11.361	17.025	1.00	41.59	O
ATOM	7970	C	ASP A 522	-63.760	14.780	17.591	1.00	37.92	C
ATOM	7971	O	ASP A 522	-63.908	15.608	18.474	1.00	38.01	O
ATOM	7973	N	GLU A 523	-64.799	14.239	16.974	1.00	37.79	N
ATOM	7974	CA	GLU A 523	-66.150	14.818	17.121	1.00	37.83	C
ATOM	7976	CB	GLU A 523	-66.618	14.919	18.591	1.00	37.85	C
ATOM	7979	CG	GLU A 523	-66.346	13.678	19.483	1.00	39.40	C
ATOM	7982	CD	GLU A 523	-66.916	13.802	20.931	1.00	41.49	C
ATOM	7983	OE1	GLU A 523	-67.714	14.737	21.207	1.00	42.81	O
ATOM	7984	OE2	GLU A 523	-66.566	12.957	21.800	1.00	42.04	O
ATOM	7985	C	GLU A 523	-66.242	16.203	16.430	1.00	37.51	C
ATOM	7986	O	GLU A 523	-67.338	16.751	16.283	1.00	37.70	O
ATOM	7988	N	LEU A 524	-65.097	16.775	16.043	1.00	37.03	N
ATOM	7989	CA	LEU A 524	-65.044	17.863	15.067	1.00	36.58	C
ATOM	7991	CB	LEU A 524	-64.174	19.022	15.559	1.00	36.56	C
ATOM	7994	CG	LEU A 524	-63.963	20.225	14.624	1.00	36.35	C
ATOM	7996	CD1	LEU A 524	-63.400	21.376	15.418	1.00	35.74	C
ATOM	8000	CD2	LEU A 524	-65.242	20.681	13.910	1.00	36.53	C
ATOM	8004	C	LEU A 524	-64.478	17.287	13.782	1.00	36.20	C
ATOM	8005	O	LEU A 524	-65.106	17.363	12.733	1.00	36.14	O
ATOM	8007	N	THR A 525	-63.298	16.686	13.869	1.00	35.83	N
ATOM	8008	CA	THR A 525	-62.756	15.935	12.748	1.00	35.66	C
ATOM	8010	CB	THR A 525	-61.585	15.015	13.173	1.00	35.48	C
ATOM	8012	OG1	THR A 525	-60.776	15.678	14.152	1.00	35.55	O
ATOM	8014	CG2	THR A 525	-60.720	14.641	11.978	1.00	34.73	C
ATOM	8018	C	THR A 525	-63.875	15.098	12.123	1.00	35.89	C
ATOM	8019	O	THR A 525	-64.040	15.095	10.909	1.00	35.98	O
ATOM	8021	N	ARG A 526	-64.667	14.420	12.954	1.00	36.06	N
ATOM	8022	CA	ARG A 526	-65.739	13.574	12.448	1.00	36.18	C
ATOM	8024	CB	ARG A 526	-66.340	12.718	13.555	1.00	36.56	C
ATOM	8027	CG	ARG A 526	-67.416	11.754	13.054	1.00	38.35	C
ATOM	8030	CD	ARG A 526	-67.781	10.702	14.092	1.00	41.00	C
ATOM	8033	NE	ARG A 526	-67.740	11.207	15.470	1.00	43.08	N
ATOM	8035	CZ	ARG A 526	-68.637	12.033	16.024	1.00	44.51	C

ATOM	8036	NH1 ARG A 526	-69.681	12.500	15.336	1.00	44.46	N
ATOM	8039	NH2 ARG A 526	-68.476	12.406	17.288	1.00	45.02	N
ATOM	8042	C ARG A 526	-66.834	14.384	11.774	1.00	35.68	C
ATOM	8043	O ARG A 526	-67.253	14.035	10.675	1.00	35.84	O
ATOM	8045	N LYS A 527	-67.309	15.447	12.424	1.00	35.06	N
ATOM	8046	CA LYS A 527	-68.285	16.343	11.780	1.00	34.61	C
ATOM	8048	CB LYS A 527	-68.680	17.528	12.683	1.00	34.74	C
ATOM	8051	CG LYS A 527	-69.818	17.225	13.651	1.00	35.40	C
ATOM	8054	CD LYS A 527	-70.301	18.452	14.441	1.00	35.99	C
ATOM	8057	CE LYS A 527	-71.280	18.021	15.556	1.00	36.42	C
ATOM	8060	NZ LYS A 527	-71.652	19.095	16.530	1.00	36.31	N
ATOM	8064	C LYS A 527	-67.723	16.860	10.457	1.00	33.77	C
ATOM	8065	O LYS A 527	-68.377	16.763	9.423	1.00	33.75	O
ATOM	8067	N ARG A 528	-66.501	17.388	10.499	1.00	32.71	N
ATOM	8068	CA ARG A 528	-65.840	17.929	9.310	1.00	31.79	C
ATOM	8070	CB ARG A 528	-64.425	18.464	9.644	1.00	31.74	C
ATOM	8073	CG ARG A 528	-64.419	19.844	10.330	1.00	30.90	C
ATOM	8076	CD ARG A 528	-63.021	20.439	10.527	1.00	29.75	C
ATOM	8079	NE ARG A 528	-63.098	21.745	11.190	1.00	29.55	N
ATOM	8081	CZ ARG A 528	-62.056	22.513	11.528	1.00	29.96	C
ATOM	8082	NH1 ARG A 528	-60.801	22.145	11.284	1.00	29.74	N
ATOM	8085	NH2 ARG A 528	-62.269	23.677	12.130	1.00	30.71	N
ATOM	8088	C ARG A 528	-65.796	16.904	8.172	1.00	30.99	C
ATOM	8089	O ARG A 528	-66.226	17.206	7.068	1.00	30.97	O
ATOM	8091	N VAL A 529	-65.313	15.694	8.450	1.00	30.12	N
ATOM	8092	CA VAL A 529	-65.239	14.631	7.438	1.00	29.43	C
ATOM	8094	CB VAL A 529	-64.557	13.354	7.972	1.00	29.39	C
ATOM	8096	CG1 VAL A 529	-64.785	12.187	7.023	1.00	28.89	C
ATOM	8100	CG2 VAL A 529	-63.062	13.589	8.186	1.00	29.44	C
ATOM	8104	C VAL A 529	-66.609	14.243	6.887	1.00	28.95	C
ATOM	8105	O VAL A 529	-66.755	14.016	5.690	1.00	29.09	O
ATOM	8107	N LEU A 530	-67.615	14.153	7.744	1.00	28.18	N
ATOM	8108	CA LEU A 530	-68.956	13.874	7.246	1.00	27.72	C
ATOM	8110	CB LEU A 530	-69.971	13.715	8.395	1.00	27.75	C
ATOM	8113	CG LEU A 530	-70.334	12.274	8.775	1.00	27.36	C
ATOM	8115	CD1 LEU A 530	-69.096	11.415	9.050	1.00	26.86	C
ATOM	8119	CD2 LEU A 530	-71.273	12.279	9.968	1.00	27.14	C
ATOM	8123	C LEU A 530	-69.401	14.963	6.256	1.00	27.16	C
ATOM	8124	O LEU A 530	-69.861	14.652	5.161	1.00	27.42	O
ATOM	8126	N SER A 531	-69.230	16.229	6.635	1.00	26.21	N
ATOM	8127	CA SER A 531	-69.730	17.359	5.853	1.00	25.12	C
ATOM	8129	CB SER A 531	-69.519	18.657	6.618	1.00	25.01	C
ATOM	8132	OG SER A 531	-68.167	18.791	6.985	1.00	23.94	O
ATOM	8134	C SER A 531	-69.055	17.476	4.507	1.00	24.48	C
ATOM	8135	O SER A 531	-69.687	17.886	3.538	1.00	24.14	O
ATOM	8137	N VAL A 532	-67.774	17.115	4.461	1.00	23.89	N
ATOM	8138	CA VAL A 532	-66.962	17.226	3.246	1.00	23.53	C
ATOM	8140	CB VAL A 532	-65.470	17.437	3.575	1.00	23.09	C
ATOM	8142	CG1 VAL A 532	-64.633	17.355	2.348	1.00	22.32	C
ATOM	8146	CG2 VAL A 532	-65.268	18.773	4.198	1.00	23.07	C
ATOM	8150	C VAL A 532	-67.100	16.028	2.316	1.00	23.67	C
ATOM	8151	O VAL A 532	-67.137	16.209	1.099	1.00	23.58	O
ATOM	8153	N ILE A 533	-67.187	14.824	2.886	1.00	23.85	N
ATOM	8154	CA ILE A 533	-67.172	13.578	2.107	1.00	24.08	C
ATOM	8156	CB ILE A 533	-66.172	12.555	2.691	1.00	24.10	C
ATOM	8158	CG1 ILE A 533	-64.745	13.079	2.631	1.00	23.47	C
ATOM	8161	CD1 ILE A 533	-64.178	13.083	1.255	1.00	23.42	C
ATOM	8165	CG2 ILE A 533	-66.254	11.231	1.932	1.00	24.56	C
ATOM	8169	C ILE A 533	-68.522	12.859	1.990	1.00	24.27	C
ATOM	8170	O ILE A 533	-68.991	12.620	.887	1.00	24.34	O

ATOM	8172	N	THR A 534	-69.133	12.487	3.111	1.00	24.57	N
ATOM	8173	CA	THR A 534	-70.279	11.564	3.072	1.00	24.92	C
ATOM	8175	CB	THR A 534	-70.200	10.522	4.207	1.00	24.93	C
ATOM	8177	OG1	THR A 534	-70.491	11.149	5.458	1.00	25.28	O
ATOM	8179	CG2	THR A 534	-68.804	9.885	4.257	1.00	24.75	C
ATOM	8183	C	THR A 534	-71.673	12.209	3.083	1.00	25.00	C
ATOM	8184	O	THR A 534	-72.601	11.646	2.522	1.00	24.94	O
ATOM	8186	N	GLU A 535	-71.821	13.374	3.708	1.00	25.26	N
ATOM	8187	CA	GLU A 535	-73.132	14.024	3.833	1.00	25.34	C
ATOM	8189	CB	GLU A 535	-73.333	14.523	5.255	1.00	25.44	C
ATOM	8192	CG	GLU A 535	-73.753	13.424	6.196	1.00	25.80	C
ATOM	8195	CD	GLU A 535	-74.307	13.963	7.468	1.00	25.72	C
ATOM	8196	OE1	GLU A 535	-75.375	13.487	7.881	1.00	24.50	O
ATOM	8197	OE2	GLU A 535	-73.678	14.877	8.039	1.00	27.03	O
ATOM	8198	C	GLU A 535	-73.363	15.179	2.859	1.00	25.26	C
ATOM	8199	O	GLU A 535	-72.686	16.206	2.939	1.00	25.02	O
ATOM	8201	N	PRO A 536	-74.351	15.030	1.959	1.00	25.25	N
ATOM	8202	CA	PRO A 536	-74.655	16.111	1.050	1.00	25.18	C
ATOM	8204	CB	PRO A 536	-75.661	15.484	.076	1.00	25.15	C
ATOM	8207	CG	PRO A 536	-75.825	14.077	.476	1.00	24.98	C
ATOM	8210	CD	PRO A 536	-75.363	13.965	1.868	1.00	25.26	C
ATOM	8213	C	PRO A 536	-75.289	17.265	1.801	1.00	25.30	C
ATOM	8214	O	PRO A 536	-75.826	17.070	2.883	1.00	25.57	O
ATOM	8215	N	ILE A 537	-75.213	18.458	1.230	1.00	25.36	N
ATOM	8216	CA	ILE A 537	-75.807	19.638	1.827	1.00	25.34	C
ATOM	8218	CB	ILE A 537	-75.201	20.918	1.221	1.00	25.14	C
ATOM	8220	CG1	ILE A 537	-73.787	21.131	1.744	1.00	24.32	C
ATOM	8223	CD1	ILE A 537	-73.228	22.472	1.406	1.00	23.22	C
ATOM	8227	CG2	ILE A 537	-76.030	22.131	1.569	1.00	25.54	C
ATOM	8231	C	ILE A 537	-77.313	19.603	1.611	1.00	25.67	C
ATOM	8232	O	ILE A 537	-77.786	19.132	.579	1.00	25.47	O
ATOM	8234	N	LEU A 538	-78.071	20.088	2.586	1.00	26.22	N
ATOM	8235	CA	LEU A 538	-79.514	20.051	2.467	1.00	26.65	C
ATOM	8237	CB	LEU A 538	-80.213	20.525	3.749	1.00	26.78	C
ATOM	8240	CG	LEU A 538	-80.162	19.609	4.989	1.00	26.64	C
ATOM	8242	CD1	LEU A 538	-81.218	20.033	6.007	1.00	26.16	C
ATOM	8246	CD2	LEU A 538	-80.346	18.138	4.633	1.00	26.12	C
ATOM	8250	C	LEU A 538	-79.925	20.881	1.262	1.00	27.08	C
ATOM	8253	N	PRO A 539	-80.786	20.306	.414	1.00	28.03	N
ATOM	8254	CA	PRO A 539	-81.117	20.855	-.888	1.00	28.44	C
ATOM	8256	CB	PRO A 539	-82.093	19.828	-1.449	1.00	28.40	C
ATOM	8259	CG	PRO A 539	-82.793	19.328	-.253	1.00	28.30	C
ATOM	8262	CD	PRO A 539	-81.720	19.228	.787	1.00	28.11	C
ATOM	8265	C	PRO A 539	-81.813	22.194	-.804	1.00	28.89	C
ATOM	8266	O	PRO A 539	-82.396	22.535	.226	1.00	28.97	O
ATOM	8267	N	PHE A 540	-81.774	22.932	-1.906	1.00	29.32	N
ATOM	8268	CA	PHE A 540	-82.380	24.245	-1.962	1.00	29.50	C
ATOM	8270	CB	PHE A 540	-82.146	24.874	-3.326	1.00	29.56	C
ATOM	8273	CG	PHE A 540	-82.757	26.227	-3.463	1.00	29.75	C
ATOM	8274	CD1	PHE A 540	-83.874	26.429	-4.255	1.00	29.47	C
ATOM	8276	CE1	PHE A 540	-84.441	27.681	-4.364	1.00	29.46	C
ATOM	8278	CZ	PHE A 540	-83.901	28.743	-3.674	1.00	29.81	C
ATOM	8280	CE2	PHE A 540	-82.791	28.553	-2.877	1.00	30.21	C
ATOM	8282	CD2	PHE A 540	-82.229	27.300	-2.770	1.00	30.17	C
ATOM	8284	C	PHE A 540	-83.865	24.129	-1.722	1.00	29.66	C
ATOM	8285	O	PHE A 540	-84.568	23.551	-2.545	1.00	29.66	O
ATOM	8287	N	GLU A 541	-84.337	24.673	-.601	1.00	29.90	N
ATOM	8288	CA	GLU A 541	-85.761	24.605	-.233	1.00	30.06	C
ATOM	8290	CB	GLU A 541	-85.951	23.724	1.017	1.00	30.15	C
ATOM	8293	CG	GLU A 541	-87.411	23.361	1.367	1.00	30.64	C

ATOM	8296	CD	GLU A 541	-88.091	24.351	2.322	1.00	31.35	C
ATOM	8297	OE1	GLU A 541	-87.468	25.370	2.680	1.00	32.67	O
ATOM	8298	OE2	GLU A 541	-89.254	24.114	2.719	1.00	30.48	O
ATOM	8299	C	GLU A 541	-86.315	26.011	-.010	1.00	29.89	C
ATOM	8300	O	GLU A 541	-86.835	26.636	-.936	1.00	29.71	O
ATOM	8302	N	LEU B 17	-69.666	-25.325	2.227	1.00	33.20	N
ATOM	8303	CA	LEU B 17	-69.356	-25.417	.755	1.00	33.49	C
ATOM	8305	CB	LEU B 17	-70.240	-26.475	.048	1.00	33.44	C
ATOM	8308	CG	LEU B 17	-70.077	-27.986	.328	1.00	33.60	C
ATOM	8310	CD1	LEU B 17	-71.285	-28.778	-.217	1.00	32.76	C
ATOM	8314	CD2	LEU B 17	-68.763	-28.553	-.230	1.00	33.21	C
ATOM	8318	C	LEU B 17	-69.513	-24.044	.055	1.00	33.54	C
ATOM	8319	O	LEU B 17	-70.550	-23.380	.195	1.00	33.77	O
ATOM	8323	N	LEU B 18	-68.481	-23.637	-.696	1.00	33.44	N
ATOM	8324	CA	LEU B 18	-68.476	-22.362	-1.454	1.00	33.07	C
ATOM	8326	CB	LEU B 18	-67.029	-21.960	-1.840	1.00	33.21	C
ATOM	8329	CG	LEU B 18	-66.065	-21.488	-.721	1.00	34.06	C
ATOM	8331	CD1	LEU B 18	-64.607	-21.351	-1.235	1.00	34.51	C
ATOM	8335	CD2	LEU B 18	-66.516	-20.164	-.061	1.00	33.43	C
ATOM	8339	C	LEU B 18	-69.379	-22.400	-2.714	1.00	32.37	C
ATOM	8340	O	LEU B 18	-69.710	-21.355	-3.274	1.00	32.29	O
ATOM	8342	N	SER B 19	-69.765	-23.597	-3.153	1.00	31.62	N
ATOM	8343	CA	SER B 19	-70.711	-23.749	-4.253	1.00	31.16	C
ATOM	8345	CB	SER B 19	-70.443	-25.067	-5.033	1.00	31.09	C
ATOM	8348	OG	SER B 19	-70.977	-26.243	-4.421	1.00	29.36	O
ATOM	8350	C	SER B 19	-72.168	-23.663	-3.759	1.00	31.45	C
ATOM	8351	O	SER B 19	-73.076	-23.432	-4.551	1.00	31.13	O
ATOM	8353	N	SER B 20	-72.386	-23.824	-2.451	1.00	31.91	N
ATOM	8354	CA	SER B 20	-73.749	-23.900	-1.879	1.00	32.34	C
ATOM	8356	CB	SER B 20	-73.706	-24.122	-.357	1.00	32.35	C
ATOM	8359	OG	SER B 20	-73.393	-25.473	-.055	1.00	32.34	O
ATOM	8361	C	SER B 20	-74.601	-22.670	-2.204	1.00	32.69	C
ATOM	8362	O	SER B 20	-74.072	-21.600	-2.487	1.00	32.89	O
ATOM	8364	N	ASP B 21	-75.921	-22.831	-2.119	1.00	33.13	N
ATOM	8365	CA	ASP B 21	-76.874	-21.901	-2.745	1.00	33.40	C
ATOM	8367	CB	ASP B 21	-78.179	-22.628	-3.045	1.00	33.64	C
ATOM	8370	CG	ASP B 21	-77.943	-23.883	-3.828	1.00	35.35	C
ATOM	8371	OD1	ASP B 21	-77.031	-23.866	-4.684	1.00	38.26	O
ATOM	8372	OD2	ASP B 21	-78.631	-24.890	-3.590	1.00	37.55	O
ATOM	8373	C	ASP B 21	-77.141	-20.649	-1.940	1.00	33.23	C
ATOM	8374	O	ASP B 21	-78.039	-20.616	-1.106	1.00	32.91	O
ATOM	8376	N	THR B 22	-76.348	-19.621	-2.229	1.00	33.46	N
ATOM	8377	CA	THR B 22	-76.443	-18.307	-1.593	1.00	33.79	C
ATOM	8379	CB	THR B 22	-76.460	-18.378	-.016	1.00	33.74	C
ATOM	8381	OG1	THR B 22	-75.426	-19.250	.460	1.00	32.97	O
ATOM	8383	CG2	THR B 22	-77.805	-18.838	.535	1.00	33.55	C
ATOM	8387	C	THR B 22	-75.217	-17.479	-2.018	1.00	34.24	C
ATOM	8388	O	THR B 22	-74.102	-17.988	-1.892	1.00	33.95	O
ATOM	8390	N	ASP B 23	-75.361	-16.248	-2.543	1.00	34.99	N
ATOM	8391	CA	ASP B 23	-76.603	-15.620	-3.103	1.00	35.61	C
ATOM	8393	CB	ASP B 23	-77.242	-16.538	-4.170	1.00	35.51	C
ATOM	8396	CG	ASP B 23	-76.196	-17.236	-5.034	1.00	36.08	C
ATOM	8397	OD1	ASP B 23	-75.084	-16.685	-5.194	1.00	35.58	O
ATOM	8398	OD2	ASP B 23	-76.473	-18.344	-5.541	1.00	37.65	O
ATOM	8399	C	ASP B 23	-77.662	-15.083	-2.097	1.00	36.10	C
ATOM	8400	O	ASP B 23	-78.644	-15.770	-1.799	1.00	36.13	O
ATOM	8402	N	GLU B 24	-77.487	-13.845	-1.616	1.00	36.74	N
ATOM	8403	CA	GLU B 24	-78.358	-13.316	-.543	1.00	37.59	C
ATOM	8405	CB	GLU B 24	-77.765	-13.666	.838	1.00	37.88	C
ATOM	8408	CG	GLU B 24	-77.624	-15.173	1.120	1.00	38.51	C

ATOM	8411	CD	GLU	B	24	-77.314	-15.510	2.588	1.00	39.48	C
ATOM	8412	OE1	GLU	B	24	-77.171	-14.581	3.434	1.00	39.60	O
ATOM	8413	OE2	GLU	B	24	-77.221	-16.725	2.884	1.00	39.51	O
ATOM	8414	C	GLU	B	24	-78.730	-11.809	-.559	1.00	38.11	C
ATOM	8415	O	GLU	B	24	-79.840	-11.453	-.972	1.00	38.20	O
ATOM	8417	N	SER	B	25	-77.825	-10.941	-.089	1.00	38.70	N
ATOM	8418	CA	SER	B	25	-78.192	-9.560	.312	1.00	39.21	C
ATOM	8420	CB	SER	B	25	-76.995	-8.783	.916	1.00	39.29	C
ATOM	8423	OG	SER	B	25	-76.267	-8.037	-.056	1.00	39.40	O
ATOM	8425	C	SER	B	25	-78.853	-8.738	-.802	1.00	39.69	C
ATOM	8426	O	SER	B	25	-78.707	-9.050	-1.995	1.00	39.75	O
ATOM	8428	N	ILE	B	26	-79.544	-7.670	-.378	1.00	40.21	N
ATOM	8429	CA	ILE	B	26	-80.554	-6.962	-1.189	1.00	40.55	C
ATOM	8431	CB	ILE	B	26	-79.925	-6.082	-2.319	1.00	40.66	C
ATOM	8433	CG1	ILE	B	26	-79.083	-4.951	-1.700	1.00	40.73	C
ATOM	8436	CD1	ILE	B	26	-78.596	-3.877	-2.698	1.00	40.68	C
ATOM	8440	CG2	ILE	B	26	-81.018	-5.475	-3.203	1.00	40.88	C
ATOM	8444	C	ILE	B	26	-81.575	-7.999	-1.714	1.00	40.72	C
ATOM	8445	O	ILE	B	26	-81.311	-8.710	-2.690	1.00	40.80	O
ATOM	8447	N	GLU	B	27	-82.729	-8.070	-1.037	1.00	40.86	N
ATOM	8448	CA	GLU	B	27	-83.694	-9.190	-1.160	1.00	40.88	C
ATOM	8450	CB	GLU	B	27	-84.782	-9.060	-.061	1.00	40.95	C
ATOM	8453	CG	GLU	B	27	-84.235	-9.257	1.379	1.00	41.28	C
ATOM	8456	CD	GLU	B	27	-85.123	-8.657	2.484	1.00	41.74	C
ATOM	8457	OE1	GLU	B	27	-85.593	-7.504	2.344	1.00	41.14	O
ATOM	8458	OE2	GLU	B	27	-85.332	-9.337	3.514	1.00	42.36	O
ATOM	8459	C	GLU	B	27	-84.290	-9.355	-2.587	1.00	40.75	C
ATOM	8460	O	GLU	B	27	-83.664	-8.936	-3.569	1.00	40.88	O
ATOM	8462	N	VAL	B	28	-85.461	-9.994	-2.712	1.00	40.40	N
ATOM	8463	CA	VAL	B	28	-86.033	-10.361	-4.027	1.00	39.97	C
ATOM	8465	CB	VAL	B	28	-86.034	-9.157	-5.053	1.00	40.03	C
ATOM	8467	CG1	VAL	B	28	-86.494	-9.600	-6.444	1.00	39.63	C
ATOM	8471	CG2	VAL	B	28	-86.900	-7.994	-4.525	1.00	39.86	C
ATOM	8475	C	VAL	B	28	-85.335	-11.600	-4.625	1.00	39.63	C
ATOM	8476	O	VAL	B	28	-85.858	-12.204	-5.567	1.00	39.63	O
ATOM	8478	N	HIS	B	29	-84.175	-11.980	-4.069	1.00	39.21	N
ATOM	8479	CA	HIS	B	29	-83.431	-13.177	-4.504	1.00	38.81	C
ATOM	8481	CB	HIS	B	29	-81.914	-12.913	-4.665	1.00	38.97	C
ATOM	8484	CG	HIS	B	29	-81.571	-11.721	-5.512	1.00	39.99	C
ATOM	8485	ND1	HIS	B	29	-81.181	-11.825	-6.834	1.00	40.79	N
ATOM	8487	CE1	HIS	B	29	-80.931	-10.615	-7.312	1.00	40.76	C
ATOM	8489	NE2	HIS	B	29	-81.134	-9.734	-6.346	1.00	40.27	N
ATOM	8491	CD2	HIS	B	29	-81.525	-10.400	-5.209	1.00	40.17	C
ATOM	8493	C	HIS	B	29	-83.611	-14.340	-3.523	1.00	38.05	C
ATOM	8494	O	HIS	B	29	-82.782	-15.259	-3.526	1.00	38.02	O
ATOM	8496	N	LYS	B	30	-84.667	-14.321	-2.694	1.00	37.06	N
ATOM	8497	CA	LYS	B	30	-85.107	-15.559	-2.022	1.00	36.33	C
ATOM	8499	CB	LYS	B	30	-86.056	-15.286	-.836	1.00	36.27	C
ATOM	8502	CG	LYS	B	30	-85.352	-14.818	.461	1.00	36.09	C
ATOM	8505	CD	LYS	B	30	-86.355	-14.599	1.612	1.00	36.03	C
ATOM	8508	CE	LYS	B	30	-85.811	-13.659	2.701	1.00	36.24	C
ATOM	8511	NZ	LYS	B	30	-86.802	-13.356	3.789	1.00	35.65	N
ATOM	8515	C	LYS	B	30	-85.710	-16.510	-3.094	1.00	35.69	C
ATOM	8516	O	LYS	B	30	-86.706	-17.205	-2.879	1.00	35.35	O
ATOM	8518	N	ASP	B	31	-85.075	-16.484	-4.268	1.00	35.13	N
ATOM	8519	CA	ASP	B	31	-85.145	-17.519	-5.282	1.00	34.64	C
ATOM	8521	CB	ASP	B	31	-84.592	-16.962	-6.623	1.00	34.60	C
ATOM	8524	CG	ASP	B	31	-84.419	-18.032	-7.725	1.00	34.24	C
ATOM	8525	OD1	ASP	B	31	-84.914	-19.165	-7.570	1.00	34.58	O
ATOM	8526	OD2	ASP	B	31	-83.776	-17.732	-8.761	1.00	31.56	O

ATOM	8527	C	ASP B	31	-84.304	-18.675	-4.721	1.00	34.31	C
ATOM	8528	O	ASP B	31	-83.122	-18.838	-5.038	1.00	34.03	O
ATOM	8530	N	LYS B	32	-84.917	-19.429	-3.813	1.00	33.94	N
ATOM	8531	CA	LYS B	32	-84.380	-20.710	-3.369	1.00	33.36	C
ATOM	8533	CB	LYS B	32	-83.772	-20.607	-1.973	1.00	33.45	C
ATOM	8536	CG	LYS B	32	-82.540	-19.693	-1.951	1.00	33.56	C
ATOM	8539	CD	LYS B	32	-82.025	-19.455	-.532	1.00	34.17	C
ATOM	8542	CE	LYS B	32	-81.671	-17.987	-.283	1.00	34.16	C
ATOM	8545	NZ	LYS B	32	-80.746	-17.437	-1.311	1.00	34.07	N
ATOM	8549	C	LYS B	32	-85.490	-21.764	-3.499	1.00	32.57	C
ATOM	8550	O	LYS B	32	-85.908	-22.412	-2.538	1.00	32.10	O
ATOM	8552	N	ALA B	33	-85.976	-21.859	-4.738	1.00	31.66	N
ATOM	8553	CA	ALA B	33	-86.576	-23.055	-5.259	1.00	30.94	C
ATOM	8555	CB	ALA B	33	-87.247	-22.773	-6.597	1.00	30.47	C
ATOM	8559	C	ALA B	33	-85.442	-24.081	-5.403	1.00	30.55	C
ATOM	8560	O	ALA B	33	-85.684	-25.230	-5.751	1.00	31.12	O
ATOM	8562	N	LYS B	34	-84.200	-23.660	-5.161	1.00	29.73	N
ATOM	8563	CA	LYS B	34	-83.107	-24.585	-4.839	1.00	29.09	C
ATOM	8565	CB	LYS B	34	-81.776	-23.817	-4.865	1.00	29.30	C
ATOM	8568	CG	LYS B	34	-80.678	-24.474	-5.726	1.00	29.76	C
ATOM	8571	CD	LYS B	34	-79.633	-23.458	-6.283	1.00	28.74	C
ATOM	8574	CE	LYS B	34	-78.417	-24.186	-6.890	1.00	27.56	C
ATOM	8577	NZ	LYS B	34	-78.745	-25.174	-7.935	1.00	26.16	N
ATOM	8581	C	LYS B	34	-83.376	-25.320	-3.472	1.00	28.37	C
ATOM	8582	O	LYS B	34	-83.254	-24.755	-2.374	1.00	27.82	O
ATOM	8584	N	LYS B	35	-83.635	-26.620	-3.604	1.00	27.53	N
ATOM	8585	CA	LYS B	35	-84.629	-27.429	-2.851	1.00	26.77	C
ATOM	8587	CB	LYS B	35	-85.729	-26.602	-2.170	1.00	26.81	C
ATOM	8590	CG	LYS B	35	-87.089	-26.520	-2.907	1.00	27.05	C
ATOM	8593	CD	LYS B	35	-88.096	-27.604	-2.452	1.00	27.41	C
ATOM	8596	CE	LYS B	35	-89.394	-27.636	-3.300	1.00	27.21	C
ATOM	8599	NZ	LYS B	35	-90.501	-26.792	-2.762	1.00	26.28	N
ATOM	8603	C	LYS B	35	-85.257	-28.357	-3.920	1.00	26.03	C
ATOM	8604	O	LYS B	35	-85.853	-29.386	-3.626	1.00	26.00	O
ATOM	8606	N	LEU B	36	-85.171	-27.909	-5.168	1.00	24.98	N
ATOM	8607	CA	LEU B	36	-85.095	-28.783	-6.319	1.00	24.02	C
ATOM	8609	CB	LEU B	36	-85.004	-27.930	-7.594	1.00	23.83	C
ATOM	8612	CG	LEU B	36	-86.208	-27.103	-8.062	1.00	22.23	C
ATOM	8614	CD1	LEU B	36	-85.762	-25.873	-8.802	1.00	19.93	C
ATOM	8618	CD2	LEU B	36	-87.092	-27.930	-8.951	1.00	22.11	C
ATOM	8622	C	LEU B	36	-83.834	-29.659	-6.185	1.00	23.58	C
ATOM	8623	O	LEU B	36	-83.840	-30.842	-6.508	1.00	23.38	O
ATOM	8625	N	GLU B	37	-82.742	-29.046	-5.739	1.00	23.18	N
ATOM	8626	CA	GLU B	37	-81.510	-29.753	-5.396	1.00	22.88	C
ATOM	8628	CB	GLU B	37	-80.466	-28.766	-4.855	1.00	22.76	C
ATOM	8631	CG	GLU B	37	-79.038	-29.297	-4.874	1.00	22.24	C
ATOM	8634	CD	GLU B	37	-78.085	-28.505	-4.009	1.00	21.06	C
ATOM	8635	OE1	GLU B	37	-78.368	-27.331	-3.737	1.00	19.65	O
ATOM	8636	OE2	GLU B	37	-77.037	-29.059	-3.615	1.00	20.62	O
ATOM	8637	C	GLU B	37	-81.748	-30.856	-4.364	1.00	22.95	C
ATOM	8638	O	GLU B	37	-81.164	-31.933	-4.458	1.00	22.99	O
ATOM	8640	N	ALA B	38	-82.591	-30.575	-3.372	1.00	23.02	N
ATOM	8641	CA	ALA B	38	-82.973	-31.565	-2.364	1.00	22.93	C
ATOM	8643	CB	ALA B	38	-83.966	-30.959	-1.374	1.00	22.73	C
ATOM	8647	C	ALA B	38	-83.580	-32.778	-3.044	1.00	23.05	C
ATOM	8648	O	ALA B	38	-83.170	-33.909	-2.813	1.00	22.52	O
ATOM	8650	N	GLU B	39	-84.546	-32.506	-3.912	1.00	23.66	N
ATOM	8651	CA	GLU B	39	-85.273	-33.534	-4.642	1.00	24.08	C
ATOM	8653	CB	GLU B	39	-86.403	-32.891	-5.453	1.00	24.28	C
ATOM	8656	CG	GLU B	39	-87.405	-33.869	-6.075	1.00	25.14	C

ATOM	8659	CD	GLU B	39	-88.773	-33.227	-6.341	1.00	26.43	C
ATOM	8660	OE1	GLU B	39	-89.232	-32.381	-5.519	1.00	25.90	O
ATOM	8661	OE2	GLU B	39	-89.387	-33.584	-7.375	1.00	26.99	O
ATOM	8662	C	GLU B	39	-84.362	-34.359	-5.545	1.00	24.18	C
ATOM	8663	O	GLU B	39	-84.483	-35.579	-5.575	1.00	24.32	O
ATOM	8665	N	VAL B	40	-83.450	-33.710	-6.269	1.00	24.43	N
ATOM	8666	CA	VAL B	40	-82.533	-34.437	-7.153	1.00	24.59	C
ATOM	8668	CB	VAL B	40	-81.734	-33.496	-8.087	1.00	24.65	C
ATOM	8670	CG1	VAL B	40	-80.611	-34.258	-8.792	1.00	24.14	C
ATOM	8674	CG2	VAL B	40	-82.662	-32.852	-9.106	1.00	23.89	C
ATOM	8678	C	VAL B	40	-81.592	-35.304	-6.328	1.00	24.99	C
ATOM	8679	O	VAL B	40	-81.265	-36.407	-6.717	1.00	24.77	O
ATOM	8681	N	ARG B	41	-81.184	-34.811	-5.171	1.00	25.79	N
ATOM	8682	CA	ARG B	41	-80.387	-35.613	-4.248	1.00	26.75	C
ATOM	8684	CB	ARG B	41	-80.032	-34.823	-2.998	1.00	27.38	C
ATOM	8687	CG	ARG B	41	-78.568	-34.777	-2.757	1.00	30.18	C
ATOM	8690	CD	ARG B	41	-78.281	-34.610	-1.273	1.00	34.43	C
ATOM	8693	NE	ARG B	41	-76.896	-34.198	-1.016	1.00	37.07	N
ATOM	8695	CZ	ARG B	41	-76.350	-33.045	-1.416	1.00	38.54	C
ATOM	8696	NH1	ARG B	41	-75.083	-32.797	-1.112	1.00	40.04	N
ATOM	8699	NH2	ARG B	41	-77.037	-32.145	-2.128	1.00	38.45	N
ATOM	8702	C	ARG B	41	-81.081	-36.860	-3.762	1.00	26.60	C
ATOM	8703	O	ARG B	41	-80.446	-37.894	-3.627	1.00	26.89	O
ATOM	8705	N	ARG B	42	-82.363	-36.733	-3.423	1.00	26.41	N
ATOM	8706	CA	ARG B	42	-83.159	-37.870	-3.002	1.00	26.09	C
ATOM	8708	CB	ARG B	42	-84.588	-37.452	-2.674	1.00	25.99	C
ATOM	8711	CG	ARG B	42	-85.426	-38.562	-2.049	1.00	25.45	C
ATOM	8714	CD	ARG B	42	-86.783	-38.078	-1.645	1.00	24.26	C
ATOM	8717	NE	ARG B	42	-87.549	-37.602	-2.794	1.00	23.73	N
ATOM	8719	CZ	ARG B	42	-88.547	-36.722	-2.717	1.00	24.51	C
ATOM	8720	NH1	ARG B	42	-88.912	-36.195	-1.545	1.00	24.99	N
ATOM	8723	NH2	ARG B	42	-89.185	-36.354	-3.816	1.00	24.46	N
ATOM	8726	C	ARG B	42	-83.167	-38.894	-4.118	1.00	26.28	C
ATOM	8727	O	ARG B	42	-82.864	-40.054	-3.888	1.00	26.39	O
ATOM	8729	N	GLU B	43	-83.476	-38.460	-5.334	1.00	26.52	N
ATOM	8730	CA	GLU B	43	-83.476	-39.370	-6.479	1.00	27.00	C
ATOM	8732	CB	GLU B	43	-83.982	-38.676	-7.753	1.00	27.50	C
ATOM	8735	CG	GLU B	43	-85.469	-38.275	-7.699	1.00	29.78	C
ATOM	8738	CD	GLU B	43	-86.342	-39.301	-6.974	1.00	32.52	C
ATOM	8739	OE1	GLU B	43	-86.388	-40.455	-7.445	1.00	35.06	O
ATOM	8740	OE2	GLU B	43	-86.966	-38.963	-5.936	1.00	33.80	O
ATOM	8741	C	GLU B	43	-82.147	-40.059	-6.769	1.00	26.37	C
ATOM	8742	O	GLU B	43	-82.162	-41.126	-7.355	1.00	26.63	O
ATOM	8744	N	ILE B	44	-81.018	-39.472	-6.373	1.00	25.85	N
ATOM	8745	CA	ILE B	44	-79.717	-40.126	-6.555	1.00	25.55	C
ATOM	8747	CB	ILE B	44	-78.547	-39.121	-6.692	1.00	25.26	C
ATOM	8749	CG1	ILE B	44	-78.706	-38.273	-7.953	1.00	24.47	C
ATOM	8752	CD1	ILE B	44	-77.690	-37.163	-8.114	1.00	23.21	C
ATOM	8756	CG2	ILE B	44	-77.207	-39.863	-6.778	1.00	25.35	C
ATOM	8760	C	ILE B	44	-79.404	-41.109	-5.428	1.00	25.84	C
ATOM	8761	O	ILE B	44	-78.925	-42.208	-5.698	1.00	25.97	O
ATOM	8763	N	ASN B	45	-79.663	-40.708	-4.177	1.00	26.25	N
ATOM	8764	CA	ASN B	45	-79.403	-41.542	-2.983	1.00	26.34	C
ATOM	8766	CB	ASN B	45	-79.284	-40.681	-1.719	1.00	26.01	C
ATOM	8769	CG	ASN B	45	-78.072	-39.789	-1.723	1.00	24.65	C
ATOM	8770	OD1	ASN B	45	-76.957	-40.226	-1.434	1.00	22.47	O
ATOM	8771	ND2	ASN B	45	-78.288	-38.516	-2.009	1.00	23.36	N
ATOM	8774	C	ASN B	45	-80.488	-42.586	-2.731	1.00	27.18	C
ATOM	8775	O	ASN B	45	-80.374	-43.390	-1.806	1.00	27.30	O
ATOM	8777	N	ASN B	46	-81.553	-42.541	-3.527	1.00	28.17	N

ATOM	8778	CA	ASN B 46	-82.619	-43.536	-3.489	1.00	29.02	C
ATOM	8780	CB	ASN B 46	-83.558	-43.289	-4.665	1.00	29.05	C
ATOM	8783	CG	ASN B 46	-84.615	-44.345	-4.806	1.00	28.82	C
ATOM	8784	OD1	ASN B 46	-84.861	-45.132	-3.895	1.00	28.14	O
ATOM	8785	ND2	ASN B 46	-85.255	-44.369	-5.968	1.00	29.67	N
ATOM	8788	C	ASN B 46	-82.076	-44.958	-3.549	1.00	30.14	C
ATOM	8789	O	ASN B 46	-81.566	-45.408	-4.581	1.00	30.30	O
ATOM	8791	N	GLU B 47	-82.218	-45.674	-2.446	1.00	31.36	N
ATOM	8792	CA	GLU B 47	-81.575	-46.968	-2.295	1.00	32.61	C
ATOM	8794	CB	GLU B 47	-81.617	-47.430	-.823	1.00	33.09	C
ATOM	8797	CG	GLU B 47	-80.974	-46.466	.218	1.00	34.19	C
ATOM	8800	CD	GLU B 47	-81.919	-45.337	.721	1.00	34.57	C
ATOM	8801	OE1	GLU B 47	-82.911	-45.012	.029	1.00	34.53	O
ATOM	8802	OE2	GLU B 47	-81.652	-44.762	1.805	1.00	34.30	O
ATOM	8803	C	GLU B 47	-82.205	-48.046	-3.172	1.00	33.18	C
ATOM	8804	O	GLU B 47	-81.604	-49.088	-3.370	1.00	33.44	O
ATOM	8806	N	LYS B 48	-83.408	-47.808	-3.692	1.00	34.11	N
ATOM	8807	CA	LYS B 48	-84.164	-48.853	-4.396	1.00	34.85	C
ATOM	8809	CB	LYS B 48	-85.379	-49.254	-3.542	1.00	35.09	C
ATOM	8812	CG	LYS B 48	-85.003	-49.723	-2.127	1.00	36.00	C
ATOM	8815	CD	LYS B 48	-86.174	-50.338	-1.350	1.00	37.13	C
ATOM	8818	CE	LYS B 48	-87.312	-49.340	-1.094	1.00	37.90	C
ATOM	8821	NZ	LYS B 48	-86.908	-48.161	-.272	1.00	38.31	N
ATOM	8825	C	LYS B 48	-84.586	-48.487	-5.835	1.00	35.13	C
ATOM	8826	O	LYS B 48	-85.544	-49.030	-6.361	1.00	34.60	O
ATOM	8828	N	ALA B 49	-83.858	-47.576	-6.474	1.00	36.09	N
ATOM	8829	CA	ALA B 49	-84.093	-47.261	-7.886	1.00	36.87	C
ATOM	8831	CB	ALA B 49	-83.550	-45.880	-8.238	1.00	36.73	C
ATOM	8835	C	ALA B 49	-83.406	-48.325	-8.714	1.00	37.56	C
ATOM	8836	O	ALA B 49	-82.315	-48.761	-8.363	1.00	37.62	O
ATOM	8838	N	GLU B 50	-84.028	-48.754	-9.807	1.00	38.54	N
ATOM	8839	CA	GLU B 50	-83.399	-49.772	-10.652	1.00	39.40	C
ATOM	8841	CB	GLU B 50	-84.380	-50.384	-11.673	1.00	39.71	C
ATOM	8844	CG	GLU B 50	-83.959	-51.808	-12.179	1.00	41.33	C
ATOM	8847	CD	GLU B 50	-83.733	-52.860	-11.039	1.00	42.88	C
ATOM	8848	OE1	GLU B 50	-84.683	-53.628	-10.718	1.00	42.94	O
ATOM	8849	OE2	GLU B 50	-82.603	-52.921	-10.472	1.00	42.96	O
ATOM	8850	C	GLU B 50	-82.171	-49.163	-11.330	1.00	39.43	C
ATOM	8851	O	GLU B 50	-82.211	-48.017	-11.791	1.00	39.66	O
ATOM	8853	N	PHE B 51	-81.084	-49.929	-11.376	1.00	39.32	N
ATOM	8854	CA	PHE B 51	-79.763	-49.360	-11.636	1.00	39.32	C
ATOM	8856	CB	PHE B 51	-78.675	-50.429	-11.476	1.00	39.62	C
ATOM	8859	CG	PHE B 51	-78.470	-50.889	-10.032	1.00	41.82	C
ATOM	8860	CD1	PHE B 51	-78.093	-49.978	-9.032	1.00	42.93	C
ATOM	8862	CE1	PHE B 51	-77.892	-50.391	-7.705	1.00	43.30	C
ATOM	8864	CZ	PHE B 51	-78.059	-51.729	-7.359	1.00	44.33	C
ATOM	8866	CE2	PHE B 51	-78.430	-52.659	-8.341	1.00	44.75	C
ATOM	8868	CD2	PHE B 51	-78.633	-52.234	-9.677	1.00	44.06	C
ATOM	8870	C	PHE B 51	-79.646	-48.625	-12.976	1.00	38.57	C
ATOM	8871	O	PHE B 51	-79.045	-47.564	-13.046	1.00	38.50	O
ATOM	8873	N	LEU B 52	-80.262	-49.155	-14.024	1.00	37.99	N
ATOM	8874	CA	LEU B 52	-80.176	-48.539	-15.351	1.00	37.42	C
ATOM	8876	CB	LEU B 52	-80.774	-49.476	-16.399	1.00	37.66	C
ATOM	8879	CG	LEU B 52	-80.502	-49.168	-17.876	1.00	38.39	C
ATOM	8881	CD1	LEU B 52	-80.171	-50.452	-18.652	1.00	38.46	C
ATOM	8885	CD2	LEU B 52	-81.703	-48.432	-18.500	1.00	39.23	C
ATOM	8889	C	LEU B 52	-80.824	-47.143	-15.431	1.00	36.69	C
ATOM	8890	O	LEU B 52	-80.383	-46.312	-16.224	1.00	36.59	O
ATOM	8892	N	THR B 53	-81.859	-46.879	-14.625	1.00	35.87	N
ATOM	8893	CA	THR B 53	-82.409	-45.510	-14.513	1.00	34.91	C

ATOM	8895	CB	THR	B	53	-83.784	-45.421	-13.805	1.00	34.72	C
ATOM	8897	OG1	THR	B	53	-84.724	-46.269	-14.455	1.00	34.04	O
ATOM	8899	CG2	THR	B	53	-84.316	-43.981	-13.855	1.00	34.81	C
ATOM	8903	C	THR	B	53	-81.454	-44.642	-13.722	1.00	34.05	C
ATOM	8904	O	THR	B	53	-81.170	-43.513	-14.121	1.00	34.16	O
ATOM	8906	N	LEU	B	54	-80.987	-45.167	-12.593	1.00	32.84	N
ATOM	8907	CA	LEU	B	54	-80.027	-44.461	-11.772	1.00	32.16	C
ATOM	8909	CB	LEU	B	54	-79.506	-45.354	-10.662	1.00	32.27	C
ATOM	8912	CG	LEU	B	54	-78.521	-44.702	-9.698	1.00	32.56	C
ATOM	8914	CD1	LEU	B	54	-79.275	-43.760	-8.769	1.00	32.38	C
ATOM	8918	CD2	LEU	B	54	-77.755	-45.789	-8.917	1.00	33.29	C
ATOM	8922	C	LEU	B	54	-78.870	-44.006	-12.634	1.00	31.50	C
ATOM	8923	O	LEU	B	54	-78.509	-42.840	-12.607	1.00	31.98	O
ATOM	8925	N	LEU	B	55	-78.297	-44.914	-13.418	1.00	30.43	N
ATOM	8926	CA	LEU	B	55	-77.220	-44.540	-14.333	1.00	29.36	C
ATOM	8928	CB	LEU	B	55	-76.765	-45.733	-15.173	1.00	29.23	C
ATOM	8931	CG	LEU	B	55	-76.157	-46.902	-14.391	1.00	28.60	C
ATOM	8933	CD1	LEU	B	55	-75.697	-48.004	-15.348	1.00	27.60	C
ATOM	8937	CD2	LEU	B	55	-75.020	-46.440	-13.469	1.00	27.20	C
ATOM	8941	C	LEU	B	55	-77.678	-43.406	-15.230	1.00	28.60	C
ATOM	8942	O	LEU	B	55	-77.063	-42.365	-15.267	1.00	28.17	O
ATOM	8944	N	GLU	B	56	-78.786	-43.594	-15.919	1.00	28.23	N
ATOM	8945	CA	GLU	B	56	-79.326	-42.530	-16.759	1.00	28.64	C
ATOM	8947	CB	GLU	B	56	-80.567	-43.031	-17.524	1.00	29.27	C
ATOM	8950	CG	GLU	B	56	-80.229	-43.783	-18.829	1.00	31.73	C
ATOM	8953	CD	GLU	B	56	-81.265	-44.860	-19.206	1.00	35.52	C
ATOM	8954	OE1	GLU	B	56	-82.474	-44.693	-18.889	1.00	36.97	O
ATOM	8955	OE2	GLU	B	56	-80.856	-45.877	-19.825	1.00	37.55	O
ATOM	8956	C	GLU	B	56	-79.632	-41.203	-15.995	1.00	27.77	C
ATOM	8957	O	GLU	B	56	-79.561	-40.103	-16.582	1.00	27.84	O
ATOM	8959	N	LEU	B	57	-79.976	-41.303	-14.710	1.00	26.29	N
ATOM	8960	CA	LEU	B	57	-80.158	-40.120	-13.895	1.00	25.06	C
ATOM	8962	CB	LEU	B	57	-80.724	-40.462	-12.514	1.00	24.83	C
ATOM	8965	CG	LEU	B	57	-80.952	-39.298	-11.544	1.00	23.72	C
ATOM	8967	CD1	LEU	B	57	-81.909	-38.278	-12.117	1.00	21.43	C
ATOM	8971	CD2	LEU	B	57	-81.477	-39.830	-10.227	1.00	22.41	C
ATOM	8975	C	LEU	B	57	-78.801	-39.459	-13.780	1.00	24.40	C
ATOM	8976	O	LEU	B	57	-78.591	-38.373	-14.306	1.00	24.70	O
ATOM	8978	N	ILE	B	58	-77.855	-40.133	-13.144	1.00	23.56	N
ATOM	8979	CA	ILE	B	58	-76.509	-39.574	-12.991	1.00	22.91	C
ATOM	8981	CB	ILE	B	58	-75.454	-40.635	-12.599	1.00	22.49	C
ATOM	8983	CG1	ILE	B	58	-75.753	-41.251	-11.235	1.00	22.09	C
ATOM	8986	CD1	ILE	B	58	-74.936	-42.464	-10.926	1.00	20.88	C
ATOM	8990	CG2	ILE	B	58	-74.103	-39.992	-12.534	1.00	22.63	C
ATOM	8994	C	ILE	B	58	-76.062	-38.927	-14.302	1.00	22.64	C
ATOM	8995	O	ILE	B	58	-75.603	-37.796	-14.315	1.00	22.38	O
ATOM	8997	N	ASP	B	59	-76.228	-39.646	-15.404	1.00	22.56	N
ATOM	8998	CA	ASP	B	59	-75.715	-39.196	-16.684	1.00	22.76	C
ATOM	9000	CB	ASP	B	59	-75.926	-40.269	-17.757	1.00	22.93	C
ATOM	9003	CG	ASP	B	59	-75.274	-39.904	-19.088	1.00	24.58	C
ATOM	9004	OD1	ASP	B	59	-74.157	-39.322	-19.081	1.00	25.70	O
ATOM	9005	OD2	ASP	B	59	-75.897	-40.186	-20.142	1.00	27.51	O
ATOM	9006	C	ASP	B	59	-76.343	-37.863	-17.104	1.00	22.40	C
ATOM	9007	O	ASP	B	59	-75.634	-36.939	-17.520	1.00	22.40	O
ATOM	9009	N	ASN	B	60	-77.662	-37.767	-16.991	1.00	21.91	N
ATOM	9010	CA	ASN	B	60	-78.349	-36.500	-17.217	1.00	21.70	C
ATOM	9012	CB	ASN	B	60	-79.867	-36.674	-17.088	1.00	22.12	C
ATOM	9015	CG	ASN	B	60	-80.477	-37.390	-18.268	1.00	22.76	C
ATOM	9016	OD1	ASN	B	60	-80.027	-37.233	-19.393	1.00	24.28	O
ATOM	9017	ND2	ASN	B	60	-81.515	-38.172	-18.017	1.00	24.36	N

ATOM	9020	C	ASN	B	60	-77.898	-35.427	-16.234	1.00	21.17	C
ATOM	9021	O	ASN	B	60	-77.700	-34.273	-16.614	1.00	20.85	O
ATOM	9023	N	VAL	B	61	-77.745	-35.805	-14.968	1.00	20.60	N
ATOM	9024	CA	VAL	B	61	-77.382	-34.841	-13.942	1.00	20.48	C
ATOM	9026	CB	VAL	B	61	-77.277	-35.490	-12.551	1.00	20.26	C
ATOM	9028	CG1	VAL	B	61	-76.585	-34.562	-11.565	1.00	19.92	C
ATOM	9032	CG2	VAL	B	61	-78.668	-35.857	-12.049	1.00	20.66	C
ATOM	9036	C	VAL	B	61	-76.068	-34.208	-14.335	1.00	20.66	C
ATOM	9037	O	VAL	B	61	-75.871	-32.998	-14.182	1.00	20.74	O
ATOM	9039	N	GLN	B	62	-75.187	-35.048	-14.875	1.00	20.92	N
ATOM	9040	CA	GLN	B	62	-73.854	-34.636	-15.285	1.00	20.69	C
ATOM	9042	CB	GLN	B	62	-72.917	-35.847	-15.404	1.00	20.77	C
ATOM	9045	CG	GLN	B	62	-72.456	-36.369	-14.035	1.00	20.91	C
ATOM	9048	CD	GLN	B	62	-71.328	-37.383	-14.112	1.00	20.75	C
ATOM	9049	OE1	GLN	B	62	-70.512	-37.498	-13.192	1.00	20.93	O
ATOM	9050	NE2	GLN	B	62	-71.291	-38.135	-15.195	1.00	20.39	N
ATOM	9053	C	GLN	B	62	-73.900	-33.834	-16.567	1.00	20.32	C
ATOM	9054	O	GLN	B	62	-73.410	-32.712	-16.569	1.00	20.79	O
ATOM	9056	N	ARG	B	63	-74.512	-34.373	-17.625	1.00	19.87	N
ATOM	9057	CA	ARG	B	63	-74.520	-33.696	-18.940	1.00	19.81	C
ATOM	9059	CB	ARG	B	63	-75.240	-34.536	-20.013	1.00	19.75	C
ATOM	9062	CG	ARG	B	63	-74.492	-35.846	-20.345	1.00	21.69	C
ATOM	9065	CD	ARG	B	63	-75.158	-36.786	-21.368	1.00	24.89	C
ATOM	9068	NE	ARG	B	63	-74.674	-36.558	-22.741	1.00	29.32	N
ATOM	9070	CZ	ARG	B	63	-75.323	-35.878	-23.702	1.00	33.71	C
ATOM	9071	NH1	ARG	B	63	-76.537	-35.345	-23.510	1.00	36.67	N
ATOM	9074	NH2	ARG	B	63	-74.761	-35.736	-24.894	1.00	34.58	N
ATOM	9077	C	ARG	B	63	-75.099	-32.276	-18.822	1.00	19.39	C
ATOM	9078	O	ARG	B	63	-74.553	-31.310	-19.395	1.00	19.63	O
ATOM	9080	N	LEU	B	64	-76.155	-32.148	-18.017	1.00	18.62	N
ATOM	9081	CA	LEU	B	64	-76.811	-30.856	-17.753	1.00	17.96	C
ATOM	9083	CB	LEU	B	64	-78.136	-31.077	-17.008	1.00	17.97	C
ATOM	9086	CG	LEU	B	64	-79.264	-31.737	-17.809	1.00	17.08	C
ATOM	9088	CD1	LEU	B	64	-80.276	-32.354	-16.879	1.00	15.15	C
ATOM	9092	CD2	LEU	B	64	-79.904	-30.734	-18.737	1.00	15.30	C
ATOM	9096	C	LEU	B	64	-75.975	-29.847	-16.966	1.00	17.39	C
ATOM	9097	O	LEU	B	64	-76.370	-28.685	-16.825	1.00	17.38	O
ATOM	9099	N	GLY	B	65	-74.848	-30.286	-16.432	1.00	16.73	N
ATOM	9100	CA	GLY	B	65	-73.917	-29.374	-15.818	1.00	16.65	C
ATOM	9103	C	GLY	B	65	-74.077	-29.249	-14.319	1.00	16.68	C
ATOM	9104	O	GLY	B	65	-73.565	-28.302	-13.718	1.00	16.79	O
ATOM	9106	N	LEU	B	66	-74.758	-30.209	-13.702	1.00	16.60	N
ATOM	9107	CA	LEU	B	66	-74.978	-30.177	-12.264	1.00	16.58	C
ATOM	9109	CB	LEU	B	66	-76.465	-30.389	-11.957	1.00	16.35	C
ATOM	9112	CG	LEU	B	66	-77.363	-29.187	-12.234	1.00	15.57	C
ATOM	9114	CD1	LEU	B	66	-78.828	-29.602	-12.163	1.00	15.26	C
ATOM	9118	CD2	LEU	B	66	-77.063	-28.051	-11.270	1.00	13.66	C
ATOM	9122	C	LEU	B	66	-74.129	-31.203	-11.513	1.00	16.94	C
ATOM	9123	O	LEU	B	66	-74.152	-31.251	-10.279	1.00	17.50	O
ATOM	9125	N	GLY	B	67	-73.373	-32.017	-12.236	1.00	17.03	N
ATOM	9126	CA	GLY	B	67	-72.541	-33.036	-11.602	1.00	17.13	C
ATOM	9129	C	GLY	B	67	-71.642	-32.557	-10.461	1.00	17.08	C
ATOM	9130	O	GLY	B	67	-71.378	-33.317	-9.522	1.00	17.33	O
ATOM	9132	N	TYR	B	68	-71.159	-31.316	-10.541	1.00	16.83	N
ATOM	9133	CA	TYR	B	68	-70.217	-30.804	-9.552	1.00	16.64	C
ATOM	9135	CB	TYR	B	68	-69.654	-29.437	-9.951	1.00	16.25	C
ATOM	9138	CG	TYR	B	68	-70.609	-28.273	-9.802	1.00	13.51	C
ATOM	9139	CD1	TYR	B	68	-70.521	-27.411	-8.738	1.00	10.49	C
ATOM	9141	CE1	TYR	B	68	-71.407	-26.342	-8.608	1.00	10.25	C
ATOM	9143	CZ	TYR	B	68	-72.382	-26.135	-9.555	1.00	9.74	C

ATOM	9144	OH	TYR	B	68	-73.253	-25.086	-9.450	1.00	7.62	O
ATOM	9146	CE2	TYR	B	68	-72.484	-26.978	-10.625	1.00	11.14	C
ATOM	9148	CD2	TYR	B	68	-71.603	-28.037	-10.748	1.00	12.56	C
ATOM	9150	C	TYR	B	68	-70.828	-30.700	-8.172	1.00	17.79	C
ATOM	9151	O	TYR	B	68	-70.107	-30.811	-7.182	1.00	18.00	O
ATOM	9153	N	ARG	B	69	-72.146	-30.485	-8.090	1.00	18.85	N
ATOM	9154	CA	ARG	B	69	-72.784	-30.269	-6.789	1.00	19.46	C
ATOM	9156	CB	ARG	B	69	-73.708	-29.047	-6.819	1.00	19.23	C
ATOM	9159	CG	ARG	B	69	-75.030	-29.219	-7.509	1.00	18.83	C
ATOM	9162	CD	ARG	B	69	-76.053	-28.192	-6.985	1.00	17.73	C
ATOM	9165	NE	ARG	B	69	-75.642	-26.830	-7.297	1.00	16.01	N
ATOM	9167	CZ	ARG	B	69	-75.330	-25.889	-6.417	1.00	14.68	C
ATOM	9168	NH1	ARG	B	69	-75.400	-26.093	-5.112	1.00	14.99	N
ATOM	9171	NH2	ARG	B	69	-74.959	-24.707	-6.861	1.00	15.04	N
ATOM	9174	C	ARG	B	69	-73.490	-31.489	-6.234	1.00	20.50	C
ATOM	9175	O	ARG	B	69	-74.084	-31.418	-5.163	1.00	20.67	O
ATOM	9177	N	PHE	B	70	-73.416	-32.607	-6.959	1.00	22.04	N
ATOM	9178	CA	PHE	B	70	-73.853	-33.917	-6.450	1.00	23.03	C
ATOM	9180	CB	PHE	B	70	-75.081	-34.428	-7.211	1.00	22.94	C
ATOM	9183	CG	PHE	B	70	-76.236	-33.503	-7.157	1.00	21.78	C
ATOM	9184	CD1	PHE	B	70	-77.024	-33.441	-6.028	1.00	20.54	C
ATOM	9186	CE1	PHE	B	70	-78.076	-32.571	-5.959	1.00	20.32	C
ATOM	9188	CZ	PHE	B	70	-78.353	-31.752	-7.025	1.00	20.73	C
ATOM	9190	CE2	PHE	B	70	-77.564	-31.796	-8.160	1.00	21.24	C
ATOM	9192	CD2	PHE	B	70	-76.511	-32.665	-8.220	1.00	21.22	C
ATOM	9194	C	PHE	B	70	-72.765	-34.962	-6.570	1.00	24.33	C
ATOM	9195	O	PHE	B	70	-73.069	-36.142	-6.614	1.00	24.56	O
ATOM	9197	N	GLU	B	71	-71.500	-34.547	-6.624	1.00	26.01	N
ATOM	9198	CA	GLU	B	71	-70.413	-35.506	-6.804	1.00	27.04	C
ATOM	9200	CB	GLU	B	71	-69.042	-34.824	-6.885	1.00	27.48	C
ATOM	9203	CG	GLU	B	71	-67.851	-35.816	-6.794	1.00	28.94	C
ATOM	9206	CD	GLU	B	71	-66.491	-35.192	-7.081	1.00	30.38	C
ATOM	9207	OE1	GLU	B	71	-66.398	-34.192	-7.835	1.00	30.52	O
ATOM	9208	OE2	GLU	B	71	-65.501	-35.730	-6.544	1.00	32.16	O
ATOM	9209	C	GLU	B	71	-70.436	-36.531	-5.680	1.00	27.46	C
ATOM	9210	O	GLU	B	71	-70.449	-37.723	-5.937	1.00	27.67	O
ATOM	9212	N	SER	B	72	-70.464	-36.087	-4.435	1.00	28.04	N
ATOM	9213	CA	SER	B	72	-70.381	-37.044	-3.348	1.00	28.86	C
ATOM	9215	CB	SER	B	72	-70.388	-36.340	-1.983	1.00	29.09	C
ATOM	9218	OG	SER	B	72	-71.649	-35.751	-1.687	1.00	30.21	O
ATOM	9220	C	SER	B	72	-71.519	-38.058	-3.481	1.00	29.30	C
ATOM	9221	O	SER	B	72	-71.304	-39.265	-3.384	1.00	29.24	O
ATOM	9223	N	ASP	B	73	-72.719	-37.560	-3.756	1.00	30.03	N
ATOM	9224	CA	ASP	B	73	-73.897	-38.417	-3.874	1.00	30.63	C
ATOM	9226	CB	ASP	B	73	-75.184	-37.598	-4.116	1.00	30.74	C
ATOM	9229	CG	ASP	B	73	-75.419	-36.509	-3.052	1.00	31.84	C
ATOM	9230	OD1	ASP	B	73	-75.562	-36.830	-1.841	1.00	31.70	O
ATOM	9231	OD2	ASP	B	73	-75.473	-35.317	-3.443	1.00	33.81	O
ATOM	9232	C	ASP	B	73	-73.722	-39.443	-4.994	1.00	30.81	C
ATOM	9233	O	ASP	B	73	-74.130	-40.590	-4.835	1.00	31.39	O
ATOM	9235	N	ILE	B	74	-73.122	-39.031	-6.112	1.00	30.96	N
ATOM	9236	CA	ILE	B	74	-72.917	-39.912	-7.282	1.00	31.04	C
ATOM	9238	CB	ILE	B	74	-72.461	-39.114	-8.524	1.00	30.88	C
ATOM	9240	CG1	ILE	B	74	-73.585	-38.213	-9.018	1.00	30.99	C
ATOM	9243	CD1	ILE	B	74	-73.105	-37.120	-9.923	1.00	31.54	C
ATOM	9247	CG2	ILE	B	74	-72.052	-40.035	-9.642	1.00	29.85	C
ATOM	9251	C	ILE	B	74	-71.885	-41.005	-6.999	1.00	31.45	C
ATOM	9252	O	ILE	B	74	-72.098	-42.165	-7.347	1.00	31.31	O
ATOM	9254	N	ARG	B	75	-70.768	-40.628	-6.380	1.00	31.93	N
ATOM	9255	CA	ARG	B	75	-69.778	-41.601	-5.947	1.00	32.38	C

ATOM	9257	CB	ARG B 75	-68.629	-40.959	-5.137	1.00	32.87	C
ATOM	9260	CG	ARG B 75	-67.310	-40.751	-5.915	1.00	35.01	C
ATOM	9263	CD	ARG B 75	-66.089	-40.489	-5.003	1.00	37.62	C
ATOM	9266	NE	ARG B 75	-65.260	-41.687	-4.783	1.00	41.03	N
ATOM	9268	CZ	ARG B 75	-64.461	-42.253	-5.700	1.00	44.29	C
ATOM	9269	NH1	ARG B 75	-64.385	-41.751	-6.931	1.00	46.22	N
ATOM	9272	NH2	ARG B 75	-63.741	-43.341	-5.404	1.00	44.30	N
ATOM	9275	C	ARG B 75	-70.478	-42.673	-5.126	1.00	32.18	C
ATOM	9276	O	ARG B 75	-70.307	-43.849	-5.398	1.00	32.19	O
ATOM	9278	N	ARG B 76	-71.280	-42.275	-4.143	1.00	32.21	N
ATOM	9279	CA	ARG B 76	-71.975	-43.256	-3.307	1.00	32.61	C
ATOM	9281	CB	ARG B 76	-72.737	-42.582	-2.162	1.00	32.75	C
ATOM	9284	CG	ARG B 76	-71.880	-42.320	-.929	1.00	33.44	C
ATOM	9287	CD	ARG B 76	-72.720	-42.046	.334	1.00	34.16	C
ATOM	9290	NE	ARG B 76	-73.797	-41.074	.124	1.00	34.53	N
ATOM	9292	CZ	ARG B 76	-73.627	-39.765	-.063	1.00	34.38	C
ATOM	9293	NH1	ARG B 76	-72.414	-39.219	-.092	1.00	34.34	N
ATOM	9296	NH2	ARG B 76	-74.689	-38.993	-.241	1.00	34.71	N
ATOM	9299	C	ARG B 76	-72.918	-44.159	-4.112	1.00	32.81	C
ATOM	9300	O	ARG B 76	-72.863	-45.379	-3.995	1.00	32.86	O
ATOM	9302	N	ALA B 77	-73.780	-43.559	-4.923	1.00	33.10	N
ATOM	9303	CA	ALA B 77	-74.655	-44.307	-5.813	1.00	33.27	C
ATOM	9305	CB	ALA B 77	-75.321	-43.366	-6.770	1.00	33.26	C
ATOM	9309	C	ALA B 77	-73.884	-45.370	-6.578	1.00	33.80	C
ATOM	9310	O	ALA B 77	-74.235	-46.542	-6.559	1.00	33.75	O
ATOM	9312	N	LEU B 78	-72.816	-44.952	-7.242	1.00	34.77	N
ATOM	9313	CA	LEU B 78	-71.987	-45.868	-8.019	1.00	35.41	C
ATOM	9315	CB	LEU B 78	-70.845	-45.109	-8.702	1.00	35.01	C
ATOM	9318	CG	LEU B 78	-71.220	-44.092	-9.782	1.00	34.37	C
ATOM	9320	CD1	LEU B 78	-69.945	-43.502	-10.381	1.00	33.93	C
ATOM	9324	CD2	LEU B 78	-72.109	-44.692	-10.873	1.00	32.94	C
ATOM	9328	C	LEU B 78	-71.409	-46.974	-7.142	1.00	36.50	C
ATOM	9329	O	LEU B 78	-71.336	-48.124	-7.552	1.00	36.87	O
ATOM	9331	N	ASP B 79	-71.005	-46.619	-5.931	1.00	37.74	N
ATOM	9332	CA	ASP B 79	-70.308	-47.545	-5.056	1.00	38.74	C
ATOM	9334	CB	ASP B 79	-69.787	-46.810	-3.819	1.00	38.99	C
ATOM	9337	CG	ASP B 79	-68.499	-47.386	-3.313	1.00	40.08	C
ATOM	9338	OD1	ASP B 79	-68.494	-48.591	-2.976	1.00	42.32	O
ATOM	9339	OD2	ASP B 79	-67.494	-46.638	-3.261	1.00	41.24	O
ATOM	9340	C	ASP B 79	-71.207	-48.705	-4.645	1.00	39.42	C
ATOM	9341	O	ASP B 79	-70.736	-49.835	-4.505	1.00	39.58	O
ATOM	9343	N	ARG B 80	-72.495	-48.424	-4.450	1.00	40.26	N
ATOM	9344	CA	ARG B 80	-73.471	-49.470	-4.143	1.00	40.92	C
ATOM	9346	CB	ARG B 80	-74.823	-48.881	-3.716	1.00	41.46	C
ATOM	9349	CG	ARG B 80	-74.871	-48.247	-2.308	1.00	42.87	C
ATOM	9352	CD	ARG B 80	-76.316	-47.826	-1.936	1.00	44.62	C
ATOM	9355	NE	ARG B 80	-76.946	-46.968	-2.949	1.00	45.87	N
ATOM	9357	CZ	ARG B 80	-76.741	-45.651	-3.089	1.00	46.82	C
ATOM	9358	NH1	ARG B 80	-75.905	-44.989	-2.282	1.00	46.27	N
ATOM	9361	NH2	ARG B 80	-77.378	-44.988	-4.059	1.00	47.22	N
ATOM	9364	C	ARG B 80	-73.667	-50.341	-5.371	1.00	40.77	C
ATOM	9365	O	ARG B 80	-73.642	-51.563	-5.275	1.00	41.06	O
ATOM	9367	N	PHE B 81	-73.863	-49.705	-6.520	1.00	40.56	N
ATOM	9368	CA	PHE B 81	-73.986	-50.418	-7.787	1.00	40.59	C
ATOM	9370	CB	PHE B 81	-73.984	-49.426	-8.956	1.00	40.73	C
ATOM	9373	CG	PHE B 81	-73.898	-50.063	-10.323	1.00	40.84	C
ATOM	9374	CD1	PHE B 81	-75.000	-50.668	-10.893	1.00	41.39	C
ATOM	9376	CE1	PHE B 81	-74.925	-51.236	-12.167	1.00	41.67	C
ATOM	9378	CZ	PHE B 81	-73.741	-51.191	-12.877	1.00	41.33	C
ATOM	9380	CE2	PHE B 81	-72.638	-50.580	-12.326	1.00	41.19	C

ATOM	9382	CD2 PHE B 81	-72.720 -50.012 -11.057	1.00 41.39	C
ATOM	9384	C PHE B 81	-72.871 -51.441 -7.953	1.00 40.47	C
ATOM	9385	O PHE B 81	-73.134 -52.583 -8.311	1.00 40.85	O
ATOM	9387	N VAL B 82	-71.633 -51.051 -7.683	1.00 40.26	N
ATOM	9388	CA VAL B 82	-70.522 -51.986 -7.821	1.00 40.24	C
ATOM	9390	CB VAL B 82	-69.173 -51.386 -7.379	1.00 40.28	C
ATOM	9392	CG1 VAL B 82	-68.164 -52.493 -7.063	1.00 39.98	C
ATOM	9396	CG2 VAL B 82	-68.644 -50.462 -8.451	1.00 40.09	C
ATOM	9400	C VAL B 82	-70.788 -53.230 -7.005	1.00 40.16	C
ATOM	9401	O VAL B 82	-70.880 -54.313 -7.559	1.00 40.13	O
ATOM	9403	N SER B 83	-70.953 -53.056 -5.698	1.00 40.26	N
ATOM	9404	CA SER B 83	-71.070 -54.174 -4.749	1.00 40.43	C
ATOM	9406	CB SER B 83	-70.797 -53.666 -3.331	1.00 40.54	C
ATOM	9409	OG SER B 83	-71.256 -52.330 -3.191	1.00 40.89	O
ATOM	9411	C SER B 83	-72.415 -54.912 -4.831	1.00 40.34	C
ATOM	9412	O SER B 83	-73.137 -55.046 -3.845	1.00 39.99	O
ATOM	9414	N SER B 84	-72.698 -55.400 -6.038	1.00 40.57	N
ATOM	9415	CA SER B 84	-73.902 -56.154 -6.411	1.00 40.63	C
ATOM	9417	CB SER B 84	-75.154 -55.679 -5.651	1.00 40.56	C
ATOM	9420	OG SER B 84	-75.292 -54.268 -5.663	1.00 39.82	O
ATOM	9422	C SER B 84	-74.104 -56.026 -7.940	1.00 40.85	C
ATOM	9423	O SER B 84	-75.104 -55.462 -8.395	1.00 41.24	O
ATOM	9425	N GLY B 85	-73.136 -56.524 -8.720	1.00 40.75	N
ATOM	9426	CA GLY B 85	-73.201 -56.506 -10.191	1.00 40.64	C
ATOM	9429	C GLY B 85	-73.225 -55.120 -10.815	1.00 40.64	C
ATOM	9430	O GLY B 85	-74.133 -54.794 -11.590	1.00 40.32	O
ATOM	9432	N THR B 93	-74.847 -57.360 -18.759	1.00 36.30	N
ATOM	9433	CA THR B 93	-75.593 -58.129 -19.774	1.00 35.99	C
ATOM	9435	CB THR B 93	-76.251 -59.379 -19.152	1.00 35.87	C
ATOM	9437	OG1 THR B 93	-76.347 -60.389 -20.158	1.00 35.64	O
ATOM	9439	CG2 THR B 93	-77.646 -59.062 -18.548	1.00 34.95	C
ATOM	9443	C THR B 93	-76.625 -57.274 -20.583	1.00 35.90	C
ATOM	9444	O THR B 93	-77.764 -57.688 -20.843	1.00 35.74	O
ATOM	9446	N SER B 94	-76.175 -56.076 -20.963	1.00 35.58	N
ATOM	9447	CA SER B 94	-76.874 -55.150 -21.866	1.00 34.97	C
ATOM	9449	CB SER B 94	-78.074 -54.472 -21.190	1.00 35.06	C
ATOM	9452	OG SER B 94	-77.702 -53.260 -20.541	1.00 34.37	O
ATOM	9454	C SER B 94	-75.831 -54.097 -22.232	1.00 34.39	C
ATOM	9455	O SER B 94	-75.212 -53.518 -21.341	1.00 34.11	O
ATOM	9457	N LEU B 95	-75.619 -53.864 -23.522	1.00 33.77	N
ATOM	9458	CA LEU B 95	-74.490 -53.045 -23.952	1.00 33.33	C
ATOM	9460	CB LEU B 95	-74.373 -53.010 -25.475	1.00 33.30	C
ATOM	9463	CG LEU B 95	-73.154 -52.246 -26.006	1.00 33.27	C
ATOM	9465	CD1 LEU B 95	-71.909 -52.455 -25.137	1.00 32.89	C
ATOM	9469	CD2 LEU B 95	-72.864 -52.649 -27.438	1.00 33.19	C
ATOM	9473	C LEU B 95	-74.551 -51.623 -23.413	1.00 33.00	C
ATOM	9474	O LEU B 95	-73.588 -51.151 -22.809	1.00 32.83	O
ATOM	9476	N HIS B 96	-75.678 -50.949 -23.633	1.00 32.64	N
ATOM	9477	CA HIS B 96	-75.853 -49.576 -23.171	1.00 32.39	C
ATOM	9479	CB HIS B 96	-77.246 -49.069 -23.527	1.00 32.64	C
ATOM	9482	CG HIS B 96	-77.528 -47.689 -23.025	1.00 34.16	C
ATOM	9483	ND1 HIS B 96	-76.565 -46.703 -22.981	1.00 36.21	N
ATOM	9485	CE1 HIS B 96	-77.096 -45.594 -22.496	1.00 36.92	C
ATOM	9487	NE2 HIS B 96	-78.371 -45.827 -22.227	1.00 37.06	N
ATOM	9489	CD2 HIS B 96	-78.664 -47.130 -22.548	1.00 35.43	C
ATOM	9491	C HIS B 96	-75.612 -49.445 -21.665	1.00 31.78	C
ATOM	9492	O HIS B 96	-74.936 -48.523 -21.215	1.00 31.95	O
ATOM	9494	N GLY B 97	-76.156 -50.372 -20.887	1.00 31.10	N
ATOM	9495	CA GLY B 97	-75.907 -50.396 -19.446	1.00 30.52	C
ATOM	9498	C GLY B 97	-74.434 -50.532 -19.085	1.00 29.95	C

ATOM	9499	O	GLY B 97	-73.936	-49.797	-18.241	1.00	30.21	O
ATOM	9501	N	THR B 98	-73.740	-51.470	-19.726	1.00	29.15	N
ATOM	9502	CA	THR B 98	-72.314	-51.690	-19.490	1.00	28.55	C
ATOM	9504	CB	THR B 98	-71.803	-52.955	-20.248	1.00	28.48	C
ATOM	9506	OG1	THR B 98	-72.678	-54.066	-19.997	1.00	28.41	O
ATOM	9508	CG2	THR B 98	-70.405	-53.331	-19.816	1.00	27.74	C
ATOM	9512	C	THR B 98	-71.492	-50.452	-19.894	1.00	28.31	C
ATOM	9513	O	THR B 98	-70.658	-49.979	-19.127	1.00	28.13	O
ATOM	9515	N	ALA B 99	-71.750	-49.910	-21.080	1.00	28.04	N
ATOM	9516	CA	ALA B 99	-70.984	-48.758	-21.578	1.00	27.88	C
ATOM	9518	CB	ALA B 99	-71.387	-48.406	-23.014	1.00	27.78	C
ATOM	9522	C	ALA B 99	-71.127	-47.538	-20.684	1.00	27.62	C
ATOM	9523	O	ALA B 99	-70.135	-46.885	-20.361	1.00	27.84	O
ATOM	9525	N	LEU B 100	-72.357	-47.233	-20.290	1.00	27.20	N
ATOM	9526	CA	LEU B 100	-72.613	-46.088	-19.422	1.00	27.08	C
ATOM	9528	CB	LEU B 100	-74.122	-45.900	-19.230	1.00	26.95	C
ATOM	9531	CG	LEU B 100	-74.560	-44.715	-18.367	1.00	26.27	C
ATOM	9533	CD1	LEU B 100	-73.753	-43.478	-18.734	1.00	26.19	C
ATOM	9537	CD2	LEU B 100	-76.052	-44.445	-18.504	1.00	25.20	C
ATOM	9541	C	LEU B 100	-71.928	-46.221	-18.046	1.00	27.23	C
ATOM	9542	O	LEU B 100	-71.368	-45.243	-17.523	1.00	27.16	O
ATOM	9544	N	SER B 101	-71.980	-47.429	-17.472	1.00	27.00	N
ATOM	9545	CA	SER B 101	-71.475	-47.682	-16.117	1.00	26.63	C
ATOM	9547	CB	SER B 101	-72.039	-48.992	-15.568	1.00	26.49	C
ATOM	9550	OG	SER B 101	-71.758	-50.063	-16.442	1.00	26.38	O
ATOM	9552	C	SER B 101	-69.956	-47.740	-16.090	1.00	26.47	C
ATOM	9553	O	SER B 101	-69.326	-47.324	-15.121	1.00	26.35	O
ATOM	9555	N	PHE B 102	-69.381	-48.288	-17.154	1.00	26.29	N
ATOM	9556	CA	PHE B 102	-67.934	-48.332	-17.322	1.00	25.91	C
ATOM	9558	CB	PHE B 102	-67.588	-49.005	-18.652	1.00	25.96	C
ATOM	9561	CG	PHE B 102	-66.133	-49.020	-18.958	1.00	25.81	C
ATOM	9562	CD1	PHE B 102	-65.360	-50.115	-18.633	1.00	25.96	C
ATOM	9564	CE1	PHE B 102	-63.999	-50.122	-18.909	1.00	26.94	C
ATOM	9566	CZ	PHE B 102	-63.403	-49.026	-19.519	1.00	26.27	C
ATOM	9568	CE2	PHE B 102	-64.174	-47.930	-19.848	1.00	26.19	C
ATOM	9570	CD2	PHE B 102	-65.531	-47.933	-19.573	1.00	26.04	C
ATOM	9572	C	PHE B 102	-67.416	-46.911	-17.315	1.00	25.56	C
ATOM	9573	O	PHE B 102	-66.486	-46.564	-16.568	1.00	25.37	O
ATOM	9575	N	ARG B 103	-68.050	-46.095	-18.152	1.00	25.20	N
ATOM	9576	CA	ARG B 103	-67.663	-44.696	-18.322	1.00	25.09	C
ATOM	9578	CB	ARG B 103	-68.510	-44.026	-19.400	1.00	25.12	C
ATOM	9581	CG	ARG B 103	-68.194	-42.561	-19.572	1.00	25.41	C
ATOM	9584	CD	ARG B 103	-68.744	-42.035	-20.889	1.00	26.68	C
ATOM	9587	NE	ARG B 103	-70.197	-41.843	-20.877	1.00	27.49	N
ATOM	9589	CZ	ARG B 103	-70.826	-40.851	-20.246	1.00	26.56	C
ATOM	9590	NH1	ARG B 103	-70.139	-39.963	-19.531	1.00	27.11	N
ATOM	9593	NH2	ARG B 103	-72.148	-40.759	-20.310	1.00	25.20	N
ATOM	9596	C	ARG B 103	-67.785	-43.898	-17.033	1.00	24.61	C
ATOM	9597	O	ARG B 103	-66.865	-43.168	-16.659	1.00	24.67	O
ATOM	9599	N	LEU B 104	-68.927	-44.013	-16.371	1.00	23.79	N
ATOM	9600	CA	LEU B 104	-69.125	-43.263	-15.148	1.00	23.23	C
ATOM	9602	CB	LEU B 104	-70.591	-43.342	-14.693	1.00	23.01	C
ATOM	9605	CG	LEU B 104	-71.607	-42.620	-15.584	1.00	21.36	C
ATOM	9607	CD1	LEU B 104	-73.002	-42.860	-15.067	1.00	19.31	C
ATOM	9611	CD2	LEU B 104	-71.310	-41.157	-15.633	1.00	19.57	C
ATOM	9615	C	LEU B 104	-68.156	-43.781	-14.071	1.00	23.15	C
ATOM	9616	O	LEU B 104	-67.445	-43.007	-13.423	1.00	23.40	O
ATOM	9618	N	LEU B 105	-68.106	-45.091	-13.894	1.00	22.79	N
ATOM	9619	CA	LEU B 105	-67.203	-45.659	-12.916	1.00	22.40	C
ATOM	9621	CB	LEU B 105	-67.303	-47.178	-12.912	1.00	22.42	C

ATOM	9624	CG	LEU B 105	-68.505	-47.709	-12.163	1.00	21.82	C
ATOM	9626	CD1	LEU B 105	-68.810	-49.120	-12.591	1.00	21.95	C
ATOM	9630	CD2	LEU B 105	-68.201	-47.648	-10.692	1.00	22.14	C
ATOM	9634	C	LEU B 105	-65.767	-45.239	-13.196	1.00	22.33	C
ATOM	9635	O	LEU B 105	-65.049	-44.823	-12.275	1.00	22.13	O
ATOM	9637	N	ARG B 106	-65.331	-45.342	-14.451	1.00	22.08	N
ATOM	9638	CA	ARG B 106	-63.953	-44.973	-14.730	1.00	22.20	C
ATOM	9640	CB	ARG B 106	-63.521	-45.266	-16.151	1.00	22.33	C
ATOM	9643	CG	ARG B 106	-62.075	-44.827	-16.329	1.00	23.60	C
ATOM	9646	CD	ARG B 106	-61.383	-45.445	-17.513	1.00	24.71	C
ATOM	9649	NE	ARG B 106	-61.078	-46.852	-17.319	1.00	24.66	N
ATOM	9651	CZ	ARG B 106	-60.418	-47.582	-18.206	1.00	25.75	C
ATOM	9652	NH1	ARG B 106	-59.995	-47.028	-19.337	1.00	26.41	N
ATOM	9655	NH2	ARG B 106	-60.182	-48.866	-17.971	1.00	26.54	N
ATOM	9658	C	ARG B 106	-63.732	-43.504	-14.438	1.00	21.75	C
ATOM	9659	O	ARG B 106	-62.801	-43.145	-13.725	1.00	21.50	O
ATOM	9661	N	GLN B 107	-64.609	-42.680	-15.003	1.00	21.69	N
ATOM	9662	CA	GLN B 107	-64.630	-41.239	-14.794	1.00	21.64	C
ATOM	9664	CB	GLN B 107	-65.964	-40.664	-15.261	1.00	21.71	C
ATOM	9667	CG	GLN B 107	-66.169	-39.178	-14.929	1.00	21.97	C
ATOM	9670	CD	GLN B 107	-67.588	-38.723	-15.172	1.00	21.54	C
ATOM	9671	OE1	GLN B 107	-68.355	-39.365	-15.906	1.00	20.08	O
ATOM	9672	NE2	GLN B 107	-67.948	-37.600	-14.556	1.00	21.45	N
ATOM	9675	C	GLN B 107	-64.454	-40.850	-13.352	1.00	21.70	C
ATOM	9676	O	GLN B 107	-63.776	-39.865	-13.078	1.00	22.14	O
ATOM	9678	N	HIS B 108	-65.091	-41.599	-12.450	1.00	21.53	N
ATOM	9679	CA	HIS B 108	-65.049	-41.326	-11.019	1.00	21.63	C
ATOM	9681	CB	HIS B 108	-66.447	-41.535	-10.429	1.00	21.47	C
ATOM	9684	CG	HIS B 108	-67.416	-40.445	-10.752	1.00	20.91	C
ATOM	9685	ND1	HIS B 108	-67.541	-39.315	-9.976	1.00	20.66	N
ATOM	9687	CE1	HIS B 108	-68.476	-38.535	-10.490	1.00	20.90	C
ATOM	9689	NE2	HIS B 108	-68.966	-39.118	-11.567	1.00	19.67	N
ATOM	9691	CD2	HIS B 108	-68.326	-40.320	-11.747	1.00	20.76	C
ATOM	9693	C	HIS B 108	-64.024	-42.209	-10.269	1.00	22.22	C
ATOM	9694	O	HIS B 108	-64.220	-42.553	-9.104	1.00	22.17	O
ATOM	9696	N	GLY B 109	-62.950	-42.613	-10.933	1.00	22.91	N
ATOM	9697	CA	GLY B 109	-61.846	-43.288	-10.245	1.00	23.85	C
ATOM	9700	C	GLY B 109	-61.986	-44.748	-9.800	1.00	24.50	C
ATOM	9701	O	GLY B 109	-61.053	-45.306	-9.189	1.00	24.12	O
ATOM	9703	N	PHE B 110	-63.128	-45.371	-10.094	1.00	25.16	N
ATOM	9704	CA	PHE B 110	-63.309	-46.791	-9.800	1.00	25.67	C
ATOM	9706	CB	PHE B 110	-64.782	-47.187	-9.891	1.00	25.80	C
ATOM	9709	CG	PHE B 110	-65.625	-46.656	-8.772	1.00	26.12	C
ATOM	9710	CD1	PHE B 110	-65.564	-47.230	-7.510	1.00	27.13	C
ATOM	9712	CE1	PHE B 110	-66.353	-46.751	-6.468	1.00	27.42	C
ATOM	9714	CZ	PHE B 110	-67.215	-45.691	-6.692	1.00	26.93	C
ATOM	9716	CE2	PHE B 110	-67.284	-45.119	-7.954	1.00	26.49	C
ATOM	9718	CD2	PHE B 110	-66.494	-45.601	-8.981	1.00	25.96	C
ATOM	9720	C	PHE B 110	-62.505	-47.615	-10.793	1.00	25.92	C
ATOM	9721	O	PHE B 110	-62.232	-47.157	-11.898	1.00	26.41	O
ATOM	9723	N	GLU B 111	-62.134	-48.832	-10.406	1.00	26.06	N
ATOM	9724	CA	GLU B 111	-61.403	-49.724	-11.304	1.00	26.17	C
ATOM	9726	CB	GLU B 111	-60.444	-50.609	-10.511	1.00	26.50	C
ATOM	9729	CG	GLU B 111	-59.372	-51.272	-11.374	1.00	28.11	C
ATOM	9732	CD	GLU B 111	-58.607	-52.377	-10.646	1.00	30.70	C
ATOM	9733	OE1	GLU B 111	-58.660	-52.436	-9.390	1.00	31.99	O
ATOM	9734	OE2	GLU B 111	-57.948	-53.190	-11.338	1.00	32.09	O
ATOM	9735	C	GLU B 111	-62.355	-50.598	-12.130	1.00	25.60	C
ATOM	9736	O	GLU B 111	-63.116	-51.388	-11.585	1.00	25.48	O
ATOM	9738	N	VAL B 112	-62.314	-50.439	-13.447	1.00	25.24	N

ATOM	9739	CA	VAL B 112	-63.026	-51.331	-14.357	1.00	24.90	C
ATOM	9741	CB	VAL B 112	-64.308	-50.712	-14.908	1.00	24.92	C
ATOM	9743	CG1	VAL B 112	-65.268	-50.437	-13.771	1.00	25.26	C
ATOM	9747	CG2	VAL B 112	-63.999	-49.450	-15.706	1.00	24.77	C
ATOM	9751	C	VAL B 112	-62.144	-51.701	-15.522	1.00	24.66	C
ATOM	9752	O	VAL B 112	-61.217	-50.964	-15.862	1.00	23.80	O
ATOM	9754	N	SER B 113	-62.467	-52.842	-16.132	1.00	24.88	N
ATOM	9755	CA	SER B 113	-61.604	-53.494	-17.109	1.00	25.21	C
ATOM	9757	CB	SER B 113	-61.286	-54.913	-16.654	1.00	24.94	C
ATOM	9760	OG	SER B 113	-60.357	-55.523	-17.528	1.00	24.31	O
ATOM	9762	C	SER B 113	-62.233	-53.532	-18.492	1.00	25.88	C
ATOM	9763	O	SER B 113	-63.446	-53.651	-18.630	1.00	25.86	O
ATOM	9765	N	GLN B 114	-61.406	-53.438	-19.526	1.00	26.92	N
ATOM	9766	CA	GLN B 114	-61.928	-53.491	-20.880	1.00	27.83	C
ATOM	9768	CB	GLN B 114	-60.846	-53.194	-21.917	1.00	27.75	C
ATOM	9771	CG	GLN B 114	-59.559	-53.952	-21.726	1.00	28.05	C
ATOM	9774	CD	GLN B 114	-58.739	-54.077	-23.007	1.00	28.39	C
ATOM	9775	OE1	GLN B 114	-59.086	-53.515	-24.051	1.00	29.05	O
ATOM	9776	NE2	GLN B 114	-57.637	-54.808	-22.925	1.00	27.62	N
ATOM	9779	C	GLN B 114	-62.615	-54.827	-21.158	1.00	28.78	C
ATOM	9780	O	GLN B 114	-63.484	-54.905	-22.022	1.00	29.10	O
ATOM	9782	N	GLU B 115	-62.248	-55.861	-20.401	1.00	29.99	N
ATOM	9783	CA	GLU B 115	-62.888	-57.181	-20.505	1.00	30.89	C
ATOM	9785	CB	GLU B 115	-62.252	-58.202	-19.549	1.00	31.10	C
ATOM	9788	CG	GLU B 115	-60.740	-58.394	-19.697	1.00	32.35	C
ATOM	9791	CD	GLU B 115	-60.320	-58.796	-21.106	1.00	34.20	C
ATOM	9792	OE1	GLU B 115	-60.965	-59.699	-21.688	1.00	34.94	O
ATOM	9793	OE2	GLU B 115	-59.348	-58.201	-21.634	1.00	35.60	O
ATOM	9794	C	GLU B 115	-64.380	-57.120	-20.224	1.00	31.26	C
ATOM	9795	O	GLU B 115	-65.125	-57.971	-20.691	1.00	31.54	O
ATOM	9797	N	ALA B 116	-64.825	-56.127	-19.465	1.00	31.91	N
ATOM	9798	CA	ALA B 116	-66.255	-55.937	-19.259	1.00	32.69	C
ATOM	9800	CB	ALA B 116	-66.521	-54.648	-18.513	1.00	32.62	C
ATOM	9804	C	ALA B 116	-66.999	-55.950	-20.597	1.00	33.33	C
ATOM	9805	O	ALA B 116	-68.156	-56.364	-20.668	1.00	33.26	O
ATOM	9807	N	PHE B 117	-66.316	-55.521	-21.655	1.00	34.16	N
ATOM	9808	CA	PHE B 117	-66.889	-55.498	-22.996	1.00	34.98	C
ATOM	9810	CB	PHE B 117	-66.310	-54.314	-23.766	1.00	35.02	C
ATOM	9813	CG	PHE B 117	-66.868	-52.997	-23.345	1.00	35.06	C
ATOM	9814	CD1	PHE B 117	-66.061	-52.046	-22.729	1.00	35.33	C
ATOM	9816	CE1	PHE B 117	-66.576	-50.825	-22.342	1.00	35.16	C
ATOM	9818	CZ	PHE B 117	-67.907	-50.544	-22.570	1.00	35.37	C
ATOM	9820	CE2	PHE B 117	-68.723	-51.488	-23.187	1.00	35.28	C
ATOM	9822	CD2	PHE B 117	-68.201	-52.701	-23.571	1.00	34.68	C
ATOM	9824	C	PHE B 117	-66.676	-56.766	-23.830	1.00	35.86	C
ATOM	9825	O	PHE B 117	-66.934	-56.755	-25.037	1.00	35.91	O
ATOM	9827	N	SER B 118	-66.215	-57.853	-23.212	1.00	36.98	N
ATOM	9828	CA	SER B 118	-65.924	-59.086	-23.960	1.00	37.76	C
ATOM	9830	CB	SER B 118	-64.836	-59.907	-23.262	1.00	37.73	C
ATOM	9833	OG	SER B 118	-65.272	-60.346	-21.987	1.00	37.77	O
ATOM	9835	C	SER B 118	-67.178	-59.939	-24.215	1.00	38.48	C
ATOM	9836	O	SER B 118	-67.207	-60.724	-25.157	1.00	38.42	O
ATOM	9838	N	GLY B 119	-68.221	-59.765	-23.407	1.00	39.51	N
ATOM	9839	CA	GLY B 119	-69.481	-60.473	-23.629	1.00	40.51	C
ATOM	9842	C	GLY B 119	-70.287	-60.034	-24.848	1.00	41.51	C
ATOM	9843	O	GLY B 119	-71.444	-60.433	-25.002	1.00	41.53	O
ATOM	9845	N	PHE B 120	-69.688	-59.210	-25.709	1.00	42.79	N
ATOM	9846	CA	PHE B 120	-70.382	-58.621	-26.864	1.00	43.81	C
ATOM	9848	CB	PHE B 120	-70.643	-57.121	-26.601	1.00	43.76	C
ATOM	9851	CG	PHE B 120	-71.367	-56.853	-25.293	1.00	43.65	C

ATOM	9852	CD1 PHE B 120	-72.762	-56.874	-25.230	1.00	43.48	C
ATOM	9854	CE1 PHE B 120	-73.437	-56.650	-24.023	1.00	43.03	C
ATOM	9856	CZ PHE B 120	-72.717	-56.413	-22.865	1.00	43.09	C
ATOM	9858	CE2 PHE B 120	-71.323	-56.397	-22.909	1.00	43.38	C
ATOM	9860	CD2 PHE B 120	-70.655	-56.619	-24.120	1.00	43.41	C
ATOM	9862	C PHE B 120	-69.607	-58.832	-28.175	1.00	44.81	C
ATOM	9863	O PHE B 120	-69.930	-58.238	-29.205	1.00	44.65	O
ATOM	9865	N LYS B 121	-68.602	-59.705	-28.126	1.00	46.19	N
ATOM	9866	CA LYS B 121	-67.757	-60.004	-29.272	1.00	47.31	C
ATOM	9868	CB LYS B 121	-66.275	-59.928	-28.870	1.00	47.47	C
ATOM	9871	CG LYS B 121	-65.743	-58.495	-28.637	1.00	48.13	C
ATOM	9874	CD LYS B 121	-64.532	-58.439	-27.675	1.00	48.89	C
ATOM	9877	CE LYS B 121	-63.261	-59.094	-28.241	1.00	49.20	C
ATOM	9880	NZ LYS B 121	-62.541	-58.232	-29.218	1.00	49.04	N
ATOM	9884	C LYS B 121	-68.096	-61.400	-29.810	1.00	48.15	C
ATOM	9885	O LYS B 121	-68.199	-62.361	-29.043	1.00	48.36	O
ATOM	9887	N ASP B 122	-68.270	-61.510	-31.126	1.00	49.04	N
ATOM	9888	CA ASP B 122	-68.611	-62.788	-31.754	1.00	49.60	C
ATOM	9890	CB ASP B 122	-69.070	-62.596	-33.217	1.00	49.52	C
ATOM	9893	CG ASP B 122	-67.989	-62.015	-34.126	1.00	49.18	C
ATOM	9894	OD1 ASP B 122	-66.792	-62.306	-33.938	1.00	48.96	O
ATOM	9895	OD2 ASP B 122	-68.351	-61.268	-35.057	1.00	48.82	O
ATOM	9896	C ASP B 122	-67.455	-63.785	-31.651	1.00	50.31	C
ATOM	9897	O ASP B 122	-66.369	-63.441	-31.171	1.00	50.33	O
ATOM	9899	N GLN B 123	-67.705	-65.014	-32.098	1.00	51.05	N
ATOM	9900	CA GLN B 123	-66.716	-66.095	-32.057	1.00	51.61	C
ATOM	9902	CB GLN B 123	-67.263	-67.330	-32.785	1.00	51.82	C
ATOM	9905	CG GLN B 123	-68.483	-67.986	-32.117	1.00	52.32	C
ATOM	9908	CD GLN B 123	-68.117	-69.137	-31.184	1.00	52.64	C
ATOM	9909	OE1 GLN B 123	-67.146	-69.062	-30.432	1.00	53.02	O
ATOM	9910	NE2 GLN B 123	-68.904	-70.207	-31.230	1.00	52.04	N
ATOM	9913	C GLN B 123	-65.358	-65.694	-32.659	1.00	51.78	C
ATOM	9914	O GLN B 123	-64.309	-66.000	-32.090	1.00	51.67	O
ATOM	9916	N ASN B 124	-65.388	-65.004	-33.799	1.00	52.05	N
ATOM	9917	CA ASN B 124	-64.166	-64.547	-34.473	1.00	52.25	C
ATOM	9919	CB ASN B 124	-64.486	-64.042	-35.884	1.00	52.25	C
ATOM	9922	CG ASN B 124	-64.911	-65.158	-36.819	1.00	51.88	C
ATOM	9923	OD1 ASN B 124	-65.918	-65.825	-36.592	1.00	51.43	O
ATOM	9924	ND2 ASN B 124	-64.144	-65.364	-37.880	1.00	51.23	N
ATOM	9927	C ASN B 124	-63.382	-63.466	-33.716	1.00	52.43	C
ATOM	9928	O ASN B 124	-62.189	-63.287	-33.959	1.00	52.33	O
ATOM	9930	N GLY B 125	-64.051	-62.747	-32.815	1.00	52.66	N
ATOM	9931	CA GLY B 125	-63.396	-61.745	-31.968	1.00	52.68	C
ATOM	9934	C GLY B 125	-63.916	-60.329	-32.135	1.00	52.63	C
ATOM	9935	O GLY B 125	-63.539	-59.447	-31.367	1.00	52.57	O
ATOM	9937	N ASN B 126	-64.782	-60.117	-33.129	1.00	52.54	N
ATOM	9938	CA ASN B 126	-65.343	-58.793	-33.443	1.00	52.41	C
ATOM	9940	CB ASN B 126	-65.575	-58.671	-34.949	1.00	52.40	C
ATOM	9943	CG ASN B 126	-64.322	-58.951	-35.751	1.00	52.64	C
ATOM	9944	OD1 ASN B 126	-63.426	-59.665	-35.297	1.00	52.77	O
ATOM	9945	ND2 ASN B 126	-64.249	-58.390	-36.952	1.00	53.00	N
ATOM	9948	C ASN B 126	-66.656	-58.524	-32.712	1.00	52.10	C
ATOM	9949	O ASN B 126	-67.253	-59.432	-32.153	1.00	52.27	O
ATOM	9951	N PHE B 127	-67.111	-57.279	-32.724	1.00	51.61	N
ATOM	9952	CA PHE B 127	-68.362	-56.936	-32.055	1.00	51.23	C
ATOM	9954	CB PHE B 127	-68.506	-55.416	-31.905	1.00	51.21	C
ATOM	9957	CG PHE B 127	-67.702	-54.848	-30.770	1.00	50.58	C
ATOM	9958	CD1 PHE B 127	-66.513	-54.192	-31.003	1.00	50.15	C
ATOM	9960	CE1 PHE B 127	-65.777	-53.686	-29.948	1.00	50.15	C
ATOM	9962	CZ PHE B 127	-66.225	-53.839	-28.647	1.00	49.71	C

ATOM	9964	CE2 PHE B 127	-67.399 -54.496 -28.405	1.00 49.50	C
ATOM	9966	CD2 PHE B 127	-68.132 -54.999 -29.459	1.00 49.99	C
ATOM	9968	C PHE B 127	-69.537 -57.505 -32.821	1.00 50.99	C
ATOM	9969	O PHE B 127	-69.468 -57.641 -34.040	1.00 50.92	O
ATOM	9971	N LEU B 128	-70.609 -57.840 -32.106	1.00 50.83	N
ATOM	9972	CA LEU B 128	-71.799 -58.422 -32.731	1.00 50.76	C
ATOM	9974	CB LEU B 128	-72.780 -58.953 -31.669	1.00 50.80	C
ATOM	9977	CG LEU B 128	-72.341 -60.098 -30.734	1.00 50.91	C
ATOM	9979	CD1 LEU B 128	-73.483 -60.532 -29.813	1.00 50.61	C
ATOM	9983	CD2 LEU B 128	-71.818 -61.303 -31.497	1.00 50.79	C
ATOM	9987	C LEU B 128	-72.500 -57.407 -33.646	1.00 50.55	C
ATOM	9988	O LEU B 128	-73.105 -56.455 -33.171	1.00 50.34	O
ATOM	9990	N GLU B 129	-72.402 -57.628 -34.957	1.00 50.50	N
ATOM	9991	CA GLU B 129	-73.038 -56.783 -35.982	1.00 50.48	C
ATOM	9993	CB GLU B 129	-73.150 -57.551 -37.310	1.00 50.68	C
ATOM	9996	CG GLU B 129	-72.262 -57.024 -38.434	1.00 51.36	C
ATOM	9999	CD GLU B 129	-72.799 -55.745 -39.060	1.00 51.88	C
ATOM	10000	OE1 GLU B 129	-72.110 -54.707 -38.961	1.00 52.72	O
ATOM	10001	OE2 GLU B 129	-73.906 -55.773 -39.643	1.00 51.58	O
ATOM	10002	C GLU B 129	-74.428 -56.264 -35.629	1.00 50.25	C
ATOM	10003	O GLU B 129	-74.737 -55.110 -35.899	1.00 50.18	O
ATOM	10005	N ASN B 130	-75.263 -57.125 -35.047	1.00 50.06	N
ATOM	10006	CA ASN B 130	-76.678 -56.806 -34.785	1.00 49.71	C
ATOM	10008	CB ASN B 130	-77.477 -58.100 -34.550	1.00 49.69	C
ATOM	10011	CG ASN B 130	-77.053 -58.837 -33.291	1.00 49.39	C
ATOM	10012	OD1 ASN B 130	-76.357 -59.848 -33.361	1.00 48.74	O
ATOM	10013	ND2 ASN B 130	-77.468 -58.329 -32.134	1.00 48.99	N
ATOM	10016	C ASN B 130	-76.934 -55.800 -33.646	1.00 49.37	C
ATOM	10017	O ASN B 130	-78.083 -55.456 -33.366	1.00 49.29	O
ATOM	10019	N LEU B 131	-75.868 -55.336 -32.996	1.00 48.98	N
ATOM	10020	CA LEU B 131	-75.957 -54.279 -31.988	1.00 48.56	C
ATOM	10022	CB LEU B 131	-74.797 -54.392 -30.991	1.00 48.44	C
ATOM	10025	CG LEU B 131	-74.759 -55.671 -30.148	1.00 48.19	C
ATOM	10027	CD1 LEU B 131	-73.382 -55.852 -29.523	1.00 47.23	C
ATOM	10031	CD2 LEU B 131	-75.861 -55.670 -29.083	1.00 47.77	C
ATOM	10035	C LEU B 131	-75.968 -52.875 -32.604	1.00 48.29	C
ATOM	10036	O LEU B 131	-76.022 -51.895 -31.874	1.00 48.44	O
ATOM	10038	N LYS B 132	-75.927 -52.771 -33.934	1.00 47.95	N
ATOM	10039	CA LYS B 132	-76.020 -51.471 -34.621	1.00 47.64	C
ATOM	10041	CB LYS B 132	-75.548 -51.590 -36.080	1.00 47.63	C
ATOM	10044	CG LYS B 132	-76.595 -52.235 -36.998	1.00 48.19	C
ATOM	10047	CD LYS B 132	-76.111 -52.483 -38.429	1.00 48.16	C
ATOM	10050	CE LYS B 132	-77.239 -53.040 -39.290	1.00 47.23	C
ATOM	10053	NZ LYS B 132	-76.725 -53.936 -40.331	1.00 47.01	N
ATOM	10057	C LYS B 132	-77.449 -50.907 -34.597	1.00 47.13	C
ATOM	10058	O LYS B 132	-77.683 -49.772 -35.005	1.00 46.87	O
ATOM	10060	N GLU B 133	-78.403 -51.716 -34.149	1.00 46.77	N
ATOM	10061	CA GLU B 133	-79.806 -51.315 -34.107	1.00 46.49	C
ATOM	10063	CB GLU B 133	-80.691 -52.458 -34.622	1.00 46.56	C
ATOM	10066	CG GLU B 133	-80.732 -52.515 -36.155	1.00 47.01	C
ATOM	10069	CD GLU B 133	-80.737 -53.926 -36.715	1.00 47.40	C
ATOM	10070	OE1 GLU B 133	-81.500 -54.775 -36.206	1.00 47.97	O
ATOM	10071	OE2 GLU B 133	-79.983 -54.177 -37.679	1.00 47.22	O
ATOM	10072	C GLU B 133	-80.260 -50.830 -32.723	1.00 45.92	C
ATOM	10073	O GLU B 133	-81.369 -50.308 -32.600	1.00 46.14	O
ATOM	10075	N ASP B 134	-79.424 -51.002 -31.691	1.00 44.98	N
ATOM	10076	CA ASP B 134	-79.611 -50.277 -30.428	1.00 44.14	C
ATOM	10078	CB ASP B 134	-79.375 -51.162 -29.197	1.00 43.95	C
ATOM	10081	CG ASP B 134	-79.646 -50.419 -27.883	1.00 43.57	C
ATOM	10082	OD1 ASP B 134	-80.014 -49.230 -27.919	1.00 42.21	O

ATOM	10083	OD2 ASP B 134	-79.488	-51.011	-26.802	1.00	43.69	O
ATOM	10084	C ASP B 134	-78.662	-49.073	-30.426	1.00	43.47	C
ATOM	10085	O ASP B 134	-77.478	-49.196	-30.104	1.00	43.55	O
ATOM	10087	N ILE B 135	-79.198	-47.908	-30.776	1.00	42.47	N
ATOM	10088	CA ILE B 135	-78.375	-46.739	-31.033	1.00	41.65	C
ATOM	10090	CB ILE B 135	-79.080	-45.783	-32.012	1.00	41.45	C
ATOM	10092	CG1 ILE B 135	-78.877	-46.313	-33.434	1.00	41.74	C
ATOM	10095	CD1 ILE B 135	-79.174	-45.331	-34.540	1.00	42.77	C
ATOM	10099	CG2 ILE B 135	-78.533	-44.391	-31.905	1.00	41.68	C
ATOM	10103	C ILE B 135	-77.901	-46.059	-29.746	1.00	41.14	C
ATOM	10104	O ILE B 135	-76.752	-45.624	-29.676	1.00	41.11	O
ATOM	10106	N LYS B 136	-78.742	-45.998	-28.714	1.00	40.53	N
ATOM	10107	CA LYS B 136	-78.281	-45.480	-27.409	1.00	40.10	C
ATOM	10109	CB LYS B 136	-79.431	-45.361	-26.375	1.00	40.40	C
ATOM	10112	CG LYS B 136	-79.910	-46.697	-25.786	1.00	42.49	C
ATOM	10115	CD LYS B 136	-81.002	-46.574	-24.698	1.00	44.79	C
ATOM	10118	CE LYS B 136	-81.599	-47.994	-24.363	1.00	46.03	C
ATOM	10121	NZ LYS B 136	-82.231	-48.149	-22.994	1.00	46.23	N
ATOM	10125	C LYS B 136	-77.104	-46.327	-26.863	1.00	38.85	C
ATOM	10126	O LYS B 136	-76.246	-45.806	-26.143	1.00	38.91	O
ATOM	10128	N ALA B 137	-77.062	-47.617	-27.219	1.00	37.27	N
ATOM	10129	CA ALA B 137	-75.935	-48.486	-26.864	1.00	35.88	C
ATOM	10131	CB ALA B 137	-76.255	-49.952	-27.116	1.00	35.69	C
ATOM	10135	C ALA B 137	-74.708	-48.085	-27.645	1.00	34.63	C
ATOM	10136	O ALA B 137	-73.647	-47.894	-27.055	1.00	34.70	O
ATOM	10138	N ILE B 138	-74.845	-47.946	-28.965	1.00	33.12	N
ATOM	10139	CA ILE B 138	-73.693	-47.597	-29.805	1.00	31.97	C
ATOM	10141	CB ILE B 138	-74.020	-47.572	-31.296	1.00	31.50	C
ATOM	10143	CG1 ILE B 138	-74.433	-48.955	-31.774	1.00	31.28	C
ATOM	10146	CD1 ILE B 138	-73.460	-50.052	-31.417	1.00	30.86	C
ATOM	10150	CG2 ILE B 138	-72.819	-47.138	-32.073	1.00	30.61	C
ATOM	10154	C ILE B 138	-73.124	-46.245	-29.413	1.00	31.57	C
ATOM	10155	O ILE B 138	-71.909	-46.090	-29.305	1.00	31.77	O
ATOM	10157	N LEU B 139	-73.997	-45.268	-29.192	1.00	30.92	N
ATOM	10158	CA LEU B 139	-73.555	-43.962	-28.713	1.00	30.29	C
ATOM	10160	CB LEU B 139	-74.732	-43.000	-28.530	1.00	30.13	C
ATOM	10163	CG LEU B 139	-74.778	-41.890	-29.574	1.00	29.96	C
ATOM	10165	CD1 LEU B 139	-73.496	-41.060	-29.508	1.00	28.26	C
ATOM	10169	CD2 LEU B 139	-76.022	-41.019	-29.397	1.00	29.48	C
ATOM	10173	C LEU B 139	-72.835	-44.136	-27.397	1.00	29.95	C
ATOM	10174	O LEU B 139	-71.709	-43.687	-27.232	1.00	30.45	O
ATOM	10176	N SER B 140	-73.492	-44.807	-26.463	1.00	29.33	N
ATOM	10177	CA SER B 140	-72.916	-45.061	-25.150	1.00	28.84	C
ATOM	10179	CB SER B 140	-73.900	-45.878	-24.298	1.00	28.93	C
ATOM	10182	OG SER B 140	-73.527	-45.870	-22.930	1.00	30.33	O
ATOM	10184	C SER B 140	-71.555	-45.770	-25.243	1.00	27.78	C
ATOM	10185	O SER B 140	-70.637	-45.454	-24.489	1.00	27.46	O
ATOM	10187	N LEU B 141	-71.432	-46.719	-26.167	1.00	26.82	N
ATOM	10188	CA LEU B 141	-70.178	-47.433	-26.348	1.00	26.34	C
ATOM	10190	CB LEU B 141	-70.366	-48.658	-27.256	1.00	26.16	C
ATOM	10193	CG LEU B 141	-69.098	-49.472	-27.584	1.00	26.09	C
ATOM	10195	CD1 LEU B 141	-68.355	-49.960	-26.330	1.00	24.64	C
ATOM	10199	CD2 LEU B 141	-69.455	-50.643	-28.481	1.00	25.94	C
ATOM	10203	C LEU B 141	-69.124	-46.476	-26.914	1.00	25.89	C
ATOM	10204	O LEU B 141	-68.025	-46.346	-26.366	1.00	25.77	O
ATOM	10206	N TYR B 142	-69.471	-45.809	-28.008	1.00	25.41	N
ATOM	10207	CA TYR B 142	-68.609	-44.804	-28.610	1.00	25.03	C
ATOM	10209	CB TYR B 142	-69.399	-43.982	-29.617	1.00	24.80	C
ATOM	10212	CG TYR B 142	-68.736	-42.695	-30.043	1.00	24.47	C
ATOM	10213	CD1 TYR B 142	-67.761	-42.686	-31.029	1.00	24.18	C

ATOM	10215	CE1 TYR B 142	-67.166	-41.516	-31.438	1.00	24.70	C
ATOM	10217	CZ TYR B 142	-67.548	-40.315	-30.871	1.00	25.80	C
ATOM	10218	OH TYR B 142	-66.940	-39.129	-31.279	1.00	27.63	O
ATOM	10220	CE2 TYR B 142	-68.522	-40.299	-29.890	1.00	25.57	C
ATOM	10222	CD2 TYR B 142	-69.110	-41.485	-29.485	1.00	24.73	C
ATOM	10224	C TYR B 142	-68.095	-43.883	-27.537	1.00	25.17	C
ATOM	10225	O TYR B 142	-66.886	-43.732	-27.360	1.00	25.40	O
ATOM	10227	N GLU B 143	-69.038	-43.293	-26.807	1.00	25.12	N
ATOM	10228	CA GLU B 143	-68.744	-42.256	-25.828	1.00	25.10	C
ATOM	10230	CB GLU B 143	-70.046	-41.726	-25.226	1.00	25.33	C
ATOM	10233	CG GLU B 143	-70.006	-40.244	-24.849	1.00	27.57	C
ATOM	10236	CD GLU B 143	-70.321	-39.313	-26.023	1.00	30.43	C
ATOM	10237	OE1 GLU B 143	-69.343	-38.764	-26.583	1.00	32.89	O
ATOM	10238	OE2 GLU B 143	-71.526	-39.130	-26.377	1.00	30.86	O
ATOM	10239	C GLU B 143	-67.794	-42.746	-24.729	1.00	24.43	C
ATOM	10240	O GLU B 143	-67.026	-41.959	-24.179	1.00	24.06	O
ATOM	10242	N ALA B 144	-67.840	-44.050	-24.449	1.00	24.01	N
ATOM	10243	CA ALA B 144	-67.033	-44.681	-23.399	1.00	23.74	C
ATOM	10245	CB ALA B 144	-67.726	-45.943	-22.883	1.00	23.48	C
ATOM	10249	C ALA B 144	-65.631	-45.025	-23.864	1.00	23.57	C
ATOM	10250	O ALA B 144	-64.720	-45.171	-23.060	1.00	23.69	O
ATOM	10252	N SER B 145	-65.455	-45.168	-25.165	1.00	23.45	N
ATOM	10253	CA SER B 145	-64.168	-45.568	-25.701	1.00	23.44	C
ATOM	10255	CB SER B 145	-64.300	-45.856	-27.214	1.00	23.63	C
ATOM	10258	OG SER B 145	-64.822	-44.748	-27.952	1.00	24.24	O
ATOM	10260	C SER B 145	-63.065	-44.520	-25.416	1.00	23.24	C
ATOM	10261	O SER B 145	-61.877	-44.848	-25.288	1.00	23.25	O
ATOM	10263	N PHE B 146	-63.454	-43.257	-25.285	1.00	22.92	N
ATOM	10264	CA PHE B 146	-62.464	-42.189	-25.170	1.00	22.35	C
ATOM	10266	CB PHE B 146	-63.090	-40.837	-25.478	1.00	22.09	C
ATOM	10269	CG PHE B 146	-63.505	-40.717	-26.911	1.00	22.40	C
ATOM	10270	CD1 PHE B 146	-62.619	-40.229	-27.868	1.00	21.98	C
ATOM	10272	CE1 PHE B 146	-62.985	-40.156	-29.201	1.00	21.80	C
ATOM	10274	CZ PHE B 146	-64.241	-40.579	-29.594	1.00	21.91	C
ATOM	10276	CE2 PHE B 146	-65.123	-41.098	-28.649	1.00	22.29	C
ATOM	10278	CD2 PHE B 146	-64.751	-41.170	-27.323	1.00	22.34	C
ATOM	10280	C PHE B 146	-61.762	-42.208	-23.837	1.00	21.91	C
ATOM	10281	O PHE B 146	-60.669	-41.683	-23.733	1.00	22.30	O
ATOM	10283	N LEU B 147	-62.361	-42.864	-22.841	1.00	21.27	N
ATOM	10284	CA LEU B 147	-61.727	-43.050	-21.535	1.00	20.49	C
ATOM	10286	CB LEU B 147	-62.791	-43.333	-20.478	1.00	20.14	C
ATOM	10289	CG LEU B 147	-63.609	-42.109	-20.102	1.00	18.87	C
ATOM	10291	CD1 LEU B 147	-64.813	-42.020	-20.959	1.00	15.94	C
ATOM	10295	CD2 LEU B 147	-63.988	-42.189	-18.642	1.00	18.62	C
ATOM	10299	C LEU B 147	-60.678	-44.163	-21.498	1.00	20.41	C
ATOM	10300	O LEU B 147	-60.243	-44.552	-20.415	1.00	20.44	O
ATOM	10302	N ALA B 148	-60.268	-44.664	-22.662	1.00	20.28	N
ATOM	10303	CA ALA B 148	-59.330	-45.784	-22.745	1.00	20.33	C
ATOM	10305	CB ALA B 148	-59.188	-46.231	-24.194	1.00	20.15	C
ATOM	10309	C ALA B 148	-57.952	-45.462	-22.170	1.00	20.54	C
ATOM	10310	O ALA B 148	-57.435	-44.366	-22.341	1.00	20.23	O
ATOM	10312	N LEU B 149	-57.366	-46.435	-21.486	1.00	21.19	N
ATOM	10313	CA LEU B 149	-55.969	-46.354	-21.067	1.00	22.02	C
ATOM	10315	CB LEU B 149	-55.768	-47.073	-19.725	1.00	21.90	C
ATOM	10318	CG LEU B 149	-56.541	-46.496	-18.523	1.00	21.50	C
ATOM	10320	CD1 LEU B 149	-55.979	-47.006	-17.223	1.00	20.96	C
ATOM	10324	CD2 LEU B 149	-56.520	-44.974	-18.505	1.00	20.95	C
ATOM	10328	C LEU B 149	-55.054	-46.930	-22.159	1.00	22.84	C
ATOM	10329	O LEU B 149	-55.507	-47.678	-23.017	1.00	22.91	O
ATOM	10331	N GLU B 150	-53.781	-46.550	-22.163	1.00	23.93	N

ATOM	10332	CA	GLU B 150	-52.858	-47.071	-23.174	1.00	25.03	C
ATOM	10334	CB	GLU B 150	-51.472	-46.431	-23.012	1.00	25.41	C
ATOM	10337	CG	GLU B 150	-50.530	-46.608	-24.213	1.00	27.44	C
ATOM	10340	CD	GLU B 150	-49.132	-45.987	-23.992	1.00	30.04	C
ATOM	10341	OE1	GLU B 150	-48.925	-45.330	-22.943	1.00	31.41	O
ATOM	10342	OE2	GLU B 150	-48.242	-46.161	-24.867	1.00	30.34	O
ATOM	10343	C	GLU B 150	-52.787	-48.600	-23.024	1.00	25.36	C
ATOM	10344	O	GLU B 150	-52.613	-49.098	-21.911	1.00	25.58	O
ATOM	10346	N	GLY B 151	-52.968	-49.337	-24.122	1.00	25.79	N
ATOM	10347	CA	GLY B 151	-52.959	-50.816	-24.086	1.00	26.11	C
ATOM	10350	C	GLY B 151	-54.323	-51.522	-24.047	1.00	26.44	C
ATOM	10351	O	GLY B 151	-54.393	-52.747	-24.135	1.00	26.91	O
ATOM	10353	N	GLU B 152	-55.407	-50.768	-23.899	1.00	26.58	N
ATOM	10354	CA	GLU B 152	-56.757	-51.319	-23.971	1.00	26.57	C
ATOM	10356	CB	GLU B 152	-57.697	-50.531	-23.065	1.00	26.74	C
ATOM	10359	CG	GLU B 152	-57.291	-50.595	-21.596	1.00	27.63	C
ATOM	10362	CD	GLU B 152	-58.271	-49.893	-20.668	1.00	28.86	C
ATOM	10363	OE1	GLU B 152	-58.239	-50.192	-19.449	1.00	28.80	O
ATOM	10364	OE2	GLU B 152	-59.066	-49.044	-21.154	1.00	29.26	O
ATOM	10365	C	GLU B 152	-57.254	-51.303	-25.415	1.00	26.53	C
ATOM	10366	O	GLU B 152	-57.937	-50.381	-25.869	1.00	26.00	O
ATOM	10368	N	ASN B 153	-56.888	-52.354	-26.127	1.00	26.77	N
ATOM	10369	CA	ASN B 153	-57.189	-52.498	-27.545	1.00	27.05	C
ATOM	10371	CB	ASN B 153	-56.345	-53.640	-28.104	1.00	27.12	C
ATOM	10374	CG	ASN B 153	-56.770	-54.992	-27.549	1.00	27.81	C
ATOM	10375	OD1	ASN B 153	-56.540	-55.307	-26.379	1.00	27.65	O
ATOM	10376	ND2	ASN B 153	-57.426	-55.782	-28.382	1.00	29.43	N
ATOM	10379	C	ASN B 153	-58.659	-52.788	-27.873	1.00	27.09	C
ATOM	10380	O	ASN B 153	-59.053	-52.658	-29.029	1.00	27.11	O
ATOM	10382	N	ILE B 154	-59.448	-53.223	-26.884	1.00	27.17	N
ATOM	10383	CA	ILE B 154	-60.865	-53.544	-27.100	1.00	27.11	C
ATOM	10385	CB	ILE B 154	-61.454	-54.369	-25.956	1.00	26.99	C
ATOM	10387	CG1	ILE B 154	-60.811	-55.750	-25.907	1.00	27.16	C
ATOM	10390	CD1	ILE B 154	-61.309	-56.622	-24.746	1.00	27.41	C
ATOM	10394	CG2	ILE B 154	-62.953	-54.519	-26.124	1.00	26.65	C
ATOM	10398	C	ILE B 154	-61.704	-52.284	-27.229	1.00	27.43	C
ATOM	10399	O	ILE B 154	-62.718	-52.279	-27.939	1.00	27.39	O
ATOM	10401	N	LEU B 155	-61.296	-51.226	-26.522	1.00	27.76	N
ATOM	10402	CA	LEU B 155	-61.964	-49.918	-26.609	1.00	27.75	C
ATOM	10404	CB	LEU B 155	-61.578	-49.024	-25.430	1.00	27.31	C
ATOM	10407	CG	LEU B 155	-61.968	-49.544	-24.052	1.00	26.64	C
ATOM	10409	CD1	LEU B 155	-61.444	-48.631	-22.981	1.00	25.82	C
ATOM	10413	CD2	LEU B 155	-63.465	-49.674	-23.944	1.00	26.45	C
ATOM	10417	C	LEU B 155	-61.620	-49.226	-27.924	1.00	28.32	C
ATOM	10418	O	LEU B 155	-62.499	-48.672	-28.583	1.00	28.72	O
ATOM	10420	N	ASP B 156	-60.348	-49.264	-28.309	1.00	28.84	N
ATOM	10421	CA	ASP B 156	-59.934	-48.737	-29.601	1.00	29.43	C
ATOM	10423	CB	ASP B 156	-58.423	-48.928	-29.823	1.00	29.92	C
ATOM	10426	CG	ASP B 156	-57.566	-47.912	-29.045	1.00	31.35	C
ATOM	10427	OD1	ASP B 156	-57.803	-46.680	-29.162	1.00	33.23	O
ATOM	10428	OD2	ASP B 156	-56.640	-48.351	-28.324	1.00	33.08	O
ATOM	10429	C	ASP B 156	-60.719	-49.411	-30.717	1.00	29.36	C
ATOM	10430	O	ASP B 156	-61.075	-48.768	-31.697	1.00	29.33	O
ATOM	10432	N	GLU B 157	-60.981	-50.705	-30.558	1.00	29.63	N
ATOM	10433	CA	GLU B 157	-61.845	-51.457	-31.479	1.00	29.98	C
ATOM	10435	CB	GLU B 157	-61.729	-52.971	-31.234	1.00	30.30	C
ATOM	10438	CG	GLU B 157	-60.664	-53.647	-32.087	1.00	31.79	C
ATOM	10441	CD	GLU B 157	-60.075	-54.901	-31.439	1.00	33.87	C
ATOM	10442	OE1	GLU B 157	-60.793	-55.560	-30.640	1.00	35.07	O
ATOM	10443	OE2	GLU B 157	-58.894	-55.223	-31.743	1.00	33.70	O

ATOM	10444	C	GLU B 157	-63.304	-51.030	-31.349	1.00	29.53	C
ATOM	10445	O	GLU B 157	-63.999	-50.878	-32.351	1.00	29.53	O
ATOM	10447	N	ALA B 158	-63.758	-50.853	-30.112	1.00	29.11	N
ATOM	10448	CA	ALA B 158	-65.104	-50.370	-29.840	1.00	28.87	C
ATOM	10450	CB	ALA B 158	-65.307	-50.195	-28.351	1.00	28.82	C
ATOM	10454	C	ALA B 158	-65.383	-49.065	-30.567	1.00	28.73	C
ATOM	10455	O	ALA B 158	-66.485	-48.857	-31.065	1.00	28.80	O
ATOM	10457	N	LYS B 159	-64.385	-48.197	-30.642	1.00	28.65	N
ATOM	10458	CA	LYS B 159	-64.553	-46.910	-31.296	1.00	28.98	C
ATOM	10460	CB	LYS B 159	-63.439	-45.950	-30.857	1.00	29.16	C
ATOM	10463	CG	LYS B 159	-63.558	-44.536	-31.426	1.00	29.82	C
ATOM	10466	CD	LYS B 159	-62.812	-43.493	-30.592	1.00	30.96	C
ATOM	10469	CE	LYS B 159	-61.295	-43.610	-30.691	1.00	31.57	C
ATOM	10472	NZ	LYS B 159	-60.630	-42.426	-30.079	1.00	31.49	N
ATOM	10476	C	LYS B 159	-64.594	-47.041	-32.826	1.00	29.08	C
ATOM	10477	O	LYS B 159	-65.385	-46.372	-33.486	1.00	28.73	O
ATOM	10479	N	VAL B 160	-63.736	-47.891	-33.386	1.00	29.53	N
ATOM	10480	CA	VAL B 160	-63.686	-48.087	-34.836	1.00	29.89	C
ATOM	10482	CB	VAL B 160	-62.466	-48.971	-35.282	1.00	29.93	C
ATOM	10484	CG1	VAL B 160	-62.576	-49.375	-36.756	1.00	29.58	C
ATOM	10488	CG2	VAL B 160	-61.152	-48.244	-35.040	1.00	29.41	C
ATOM	10492	C	VAL B 160	-65.001	-48.723	-35.263	1.00	30.32	C
ATOM	10493	O	VAL B 160	-65.507	-48.459	-36.362	1.00	30.41	O
ATOM	10495	N	PHE B 161	-65.558	-49.540	-34.371	1.00	30.83	N
ATOM	10496	CA	PHE B 161	-66.854	-50.178	-34.599	1.00	31.34	C
ATOM	10498	CB	PHE B 161	-67.090	-51.321	-33.599	1.00	31.34	C
ATOM	10501	CG	PHE B 161	-68.492	-51.834	-33.603	1.00	31.09	C
ATOM	10502	CD1	PHE B 161	-68.940	-52.636	-34.633	1.00	31.61	C
ATOM	10504	CE1	PHE B 161	-70.249	-53.098	-34.649	1.00	31.48	C
ATOM	10506	CZ	PHE B 161	-71.114	-52.748	-33.628	1.00	30.76	C
ATOM	10508	CE2	PHE B 161	-70.677	-51.944	-32.606	1.00	30.19	C
ATOM	10510	CD2	PHE B 161	-69.378	-51.486	-32.597	1.00	30.70	C
ATOM	10512	C	PHE B 161	-67.992	-49.173	-34.504	1.00	31.65	C
ATOM	10513	O	PHE B 161	-68.785	-49.038	-35.432	1.00	31.62	O
ATOM	10515	N	ALA B 162	-68.068	-48.483	-33.373	1.00	32.26	N
ATOM	10516	CA	ALA B 162	-69.135	-47.519	-33.129	1.00	32.84	C
ATOM	10518	CB	ALA B 162	-68.948	-46.854	-31.778	1.00	32.66	C
ATOM	10522	C	ALA B 162	-69.224	-46.474	-34.245	1.00	33.40	C
ATOM	10523	O	ALA B 162	-70.229	-46.408	-34.937	1.00	33.45	O
ATOM	10525	N	ILE B 163	-68.164	-45.697	-34.444	1.00	34.38	N
ATOM	10526	CA	ILE B 163	-68.166	-44.610	-35.439	1.00	35.23	C
ATOM	10528	CB	ILE B 163	-66.734	-44.062	-35.733	1.00	35.18	C
ATOM	10530	CG1	ILE B 163	-66.092	-43.457	-34.488	1.00	35.21	C
ATOM	10533	CD1	ILE B 163	-64.620	-43.171	-34.663	1.00	35.96	C
ATOM	10537	CG2	ILE B 163	-66.778	-42.977	-36.799	1.00	34.72	C
ATOM	10541	C	ILE B 163	-68.778	-45.052	-36.766	1.00	36.18	C
ATOM	10542	O	ILE B 163	-69.588	-44.335	-37.360	1.00	35.87	O
ATOM	10544	N	SER B 164	-68.379	-46.240	-37.217	1.00	37.67	N
ATOM	10545	CA	SER B 164	-68.761	-46.752	-38.539	1.00	38.74	C
ATOM	10547	CB	SER B 164	-68.104	-48.117	-38.815	1.00	38.80	C
ATOM	10550	OG	SER B 164	-68.733	-49.158	-38.077	1.00	39.07	O
ATOM	10552	C	SER B 164	-70.277	-46.861	-38.699	1.00	39.51	C
ATOM	10553	O	SER B 164	-70.820	-46.494	-39.737	1.00	39.96	O
ATOM	10555	N	HIS B 165	-70.962	-47.363	-37.680	1.00	40.28	N
ATOM	10556	CA	HIS B 165	-72.410	-47.485	-37.767	1.00	41.07	C
ATOM	10558	CB	HIS B 165	-72.911	-48.686	-36.957	1.00	41.38	C
ATOM	10561	CG	HIS B 165	-72.571	-50.005	-37.587	1.00	42.65	C
ATOM	10562	ND1	HIS B 165	-71.647	-50.875	-37.046	1.00	43.54	N
ATOM	10564	CE1	HIS B 165	-71.535	-51.938	-37.825	1.00	43.54	C
ATOM	10566	NE2	HIS B 165	-72.346	-51.785	-38.858	1.00	43.25	N

ATOM	10568	CD2 HIS B 165	-73.001 -50.582 -38.737	1.00 43.31	C
ATOM	10570	C HIS B 165	-73.125 -46.193 -37.380	1.00 41.21	C
ATOM	10571	O HIS B 165	-74.257 -45.980 -37.805	1.00 41.53	O
ATOM	10573	N LEU B 166	-72.459 -45.326 -36.614	1.00 41.32	N
ATOM	10574	CA LEU B 166	-73.042 -44.039 -36.201	1.00 41.31	C
ATOM	10576	CB LEU B 166	-72.317 -43.460 -34.973	1.00 41.14	C
ATOM	10579	CG LEU B 166	-72.732 -43.975 -33.589	1.00 40.03	C
ATOM	10581	CD1 LEU B 166	-71.673 -43.634 -32.579	1.00 38.86	C
ATOM	10585	CD2 LEU B 166	-74.077 -43.421 -33.151	1.00 38.52	C
ATOM	10589	C LEU B 166	-73.045 -42.992 -37.313	1.00 41.80	C
ATOM	10590	O LEU B 166	-74.006 -42.230 -37.442	1.00 41.73	O
ATOM	10592	N LYS B 167	-71.982 -42.947 -38.116	1.00 42.47	N
ATOM	10593	CA LYS B 167	-71.847 -41.891 -39.137	1.00 43.00	C
ATOM	10595	CB LYS B 167	-70.416 -41.822 -39.702	1.00 43.12	C
ATOM	10598	CD LYS B 167	-69.937 -43.068 -40.450	1.00 44.12	C
ATOM	10601	CD LYS B 167	-69.138 -42.725 -41.736	1.00 45.36	C
ATOM	10604	CE LYS B 167	-67.898 -41.842 -41.483	1.00 45.65	C
ATOM	10607	NZ LYS B 167	-67.267 -41.400 -42.762	1.00 44.94	N
ATOM	10611	C LYS B 167	-72.873 -41.950 -40.283	1.00 42.99	C
ATOM	10612	O LYS B 167	-72.826 -41.124 -41.182	1.00 42.95	O
ATOM	10614	N GLU B 168	-73.791 -42.912 -40.244	1.00 43.23	N
ATOM	10615	CA GLU B 168	-74.928 -42.940 -41.167	1.00 43.44	C
ATOM	10617	CB GLU B 168	-74.652 -43.969 -42.272	1.00 43.58	C
ATOM	10620	CG GLU B 168	-73.918 -43.363 -43.501	1.00 44.36	C
ATOM	10623	CD GLU B 168	-72.618 -44.082 -43.895	1.00 44.53	C
ATOM	10624	OE1 GLU B 168	-71.853 -44.494 -42.992	1.00 44.71	O
ATOM	10625	OE2 GLU B 168	-72.351 -44.197 -45.115	1.00 43.31	O
ATOM	10626	C GLU B 168	-76.263 -43.214 -40.442	1.00 43.24	C
ATOM	10627	O GLU B 168	-76.932 -42.291 -39.942	1.00 42.54	O
ATOM	10629	N GLY B 175	-83.548 -42.037 -36.239	1.00 49.43	N
ATOM	10630	CA GLY B 175	-84.473 -41.266 -35.416	1.00 49.54	C
ATOM	10633	C GLY B 175	-84.709 -39.884 -36.001	1.00 49.76	C
ATOM	10634	O GLY B 175	-84.828 -39.741 -37.220	1.00 49.70	O
ATOM	10636	N LYS B 176	-84.766 -38.870 -35.131	1.00 49.95	N
ATOM	10637	CA LYS B 176	-85.014 -37.468 -35.534	1.00 50.05	C
ATOM	10639	CB LYS B 176	-86.522 -37.225 -35.730	1.00 50.27	C
ATOM	10642	CG LYS B 176	-87.301 -36.813 -34.460	1.00 51.33	C
ATOM	10645	CD LYS B 176	-88.553 -37.658 -34.202	1.00 52.21	C
ATOM	10648	CE LYS B 176	-88.841 -37.754 -32.695	1.00 52.36	C
ATOM	10651	NZ LYS B 176	-90.029 -38.607 -32.422	1.00 52.15	N
ATOM	10655	C LYS B 176	-84.427 -36.432 -34.554	1.00 49.75	C
ATOM	10656	O LYS B 176	-83.990 -35.370 -34.972	1.00 49.88	O
ATOM	10658	N GLU B 177	-84.473 -36.729 -33.254	1.00 49.46	N
ATOM	10659	CA GLU B 177	-83.769 -35.959 -32.221	1.00 48.92	C
ATOM	10661	CB GLU B 177	-84.628 -35.813 -30.946	1.00 49.03	C
ATOM	10664	CG GLU B 177	-84.196 -36.691 -29.732	1.00 49.60	C
ATOM	10667	CD GLU B 177	-85.278 -36.871 -28.662	1.00 50.10	C
ATOM	10668	OE1 GLU B 177	-84.998 -37.576 -27.670	1.00 49.78	O
ATOM	10669	OE2 GLU B 177	-86.401 -36.336 -28.808	1.00 50.78	O
ATOM	10670	C GLU B 177	-82.469 -36.703 -31.924	1.00 48.19	C
ATOM	10671	O GLU B 177	-81.421 -36.085 -31.721	1.00 48.68	O
ATOM	10673	N LEU B 178	-82.557 -38.037 -31.900	1.00 46.98	N
ATOM	10674	CA LEU B 178	-81.398 -38.924 -31.817	1.00 45.91	C
ATOM	10676	CB LEU B 178	-81.838 -40.387 -31.945	1.00 45.70	C
ATOM	10679	CG LEU B 178	-81.111 -41.408 -31.074	1.00 45.24	C
ATOM	10681	CD1 LEU B 178	-81.498 -41.225 -29.609	1.00 44.88	C
ATOM	10685	CD2 LEU B 178	-81.413 -42.826 -31.539	1.00 44.19	C
ATOM	10689	C LEU B 178	-80.423 -38.571 -32.934	1.00 45.27	C
ATOM	10690	O LEU B 178	-79.211 -38.624 -32.759	1.00 44.99	O
ATOM	10692	N ALA B 179	-80.970 -38.203 -34.088	1.00 44.75	N

ATOM	10693	CA	ALA B 179	-80.174	-37.661	-35.181	1.00	44.24	C
ATOM	10695	CB	ALA B 179	-81.088	-37.069	-36.231	1.00	44.30	C
ATOM	10699	C	ALA B 179	-79.192	-36.602	-34.683	1.00	43.65	C
ATOM	10700	O	ALA B 179	-78.028	-36.604	-35.061	1.00	43.56	O
ATOM	10702	N	GLU B 180	-79.675	-35.707	-33.827	1.00	43.01	N
ATOM	10703	CA	GLU B 180	-78.863	-34.610	-33.311	1.00	42.49	C
ATOM	10705	CB	GLU B 180	-79.749	-33.477	-32.789	1.00	42.85	C
ATOM	10708	CG	GLU B 180	-80.574	-32.800	-33.881	1.00	44.23	C
ATOM	10711	CD	GLU B 180	-80.991	-31.382	-33.519	1.00	45.80	C
ATOM	10712	OE1	GLU B 180	-80.098	-30.580	-33.156	1.00	46.69	O
ATOM	10713	OE2	GLU B 180	-82.203	-31.065	-33.610	1.00	46.71	O
ATOM	10714	C	GLU B 180	-77.944	-35.075	-32.210	1.00	41.39	C
ATOM	10715	O	GLU B 180	-76.843	-34.562	-32.071	1.00	41.27	O
ATOM	10717	N	GLN B 181	-78.404	-36.035	-31.419	1.00	40.20	N
ATOM	10718	CA	GLN B 181	-77.570	-36.632	-30.394	1.00	39.45	C
ATOM	10720	CB	GLN B 181	-78.280	-37.809	-29.736	1.00	39.88	C
ATOM	10723	CG	GLN B 181	-78.005	-37.959	-28.249	1.00	41.64	C
ATOM	10726	CD	GLN B 181	-78.873	-37.036	-27.404	1.00	44.00	C
ATOM	10727	OE1	GLN B 181	-80.109	-36.992	-27.581	1.00	45.65	O
ATOM	10728	NE2	GLN B 181	-78.235	-36.291	-26.474	1.00	43.41	N
ATOM	10731	C	GLN B 181	-76.270	-37.120	-31.003	1.00	38.28	C
ATOM	10732	O	GLN B 181	-75.203	-36.767	-30.523	1.00	38.27	O
ATOM	10734	N	VAL B 182	-76.358	-37.911	-32.072	1.00	37.05	N
ATOM	10735	CA	VAL B 182	-75.160	-38.515	-32.675	1.00	36.20	C
ATOM	10737	CB	VAL B 182	-75.486	-39.661	-33.681	1.00	36.11	C
ATOM	10739	CG1	VAL B 182	-76.493	-40.641	-33.080	1.00	36.13	C
ATOM	10743	CG2	VAL B 182	-75.977	-39.106	-35.004	1.00	36.06	C
ATOM	10747	C	VAL B 182	-74.238	-37.495	-33.364	1.00	35.44	C
ATOM	10748	O	VAL B 182	-73.012	-37.682	-33.375	1.00	35.19	O
ATOM	10750	N	SER B 183	-74.811	-36.432	-33.936	1.00	34.41	N
ATOM	10751	CA	SER B 183	-73.998	-35.416	-34.630	1.00	33.74	C
ATOM	10753	CB	SER B 183	-74.858	-34.438	-35.422	1.00	33.58	C
ATOM	10756	OG	SER B 183	-76.003	-35.085	-35.925	1.00	34.52	O
ATOM	10758	C	SER B 183	-73.184	-34.644	-33.625	1.00	32.88	C
ATOM	10759	O	SER B 183	-72.055	-34.251	-33.906	1.00	33.25	O
ATOM	10761	N	HIS B 184	-73.788	-34.427	-32.461	1.00	31.83	N
ATOM	10762	CA	HIS B 184	-73.150	-33.787	-31.326	1.00	30.98	C
ATOM	10764	CB	HIS B 184	-74.188	-33.627	-30.212	1.00	31.14	C
ATOM	10767	CG	HIS B 184	-73.710	-32.849	-29.030	1.00	31.91	C
ATOM	10768	ND1	HIS B 184	-73.326	-31.527	-29.114	1.00	33.18	N
ATOM	10770	CE1	HIS B 184	-72.960	-31.107	-27.914	1.00	33.10	C
ATOM	10772	NE2	HIS B 184	-73.106	-32.102	-27.056	1.00	31.85	N
ATOM	10774	CD2	HIS B 184	-73.583	-33.200	-27.727	1.00	31.86	C
ATOM	10776	C	HIS B 184	-71.968	-34.632	-30.865	1.00	30.13	C
ATOM	10777	O	HIS B 184	-70.863	-34.122	-30.709	1.00	29.95	O
ATOM	10779	N	ALA B 185	-72.194	-35.931	-30.683	1.00	29.27	N
ATOM	10780	CA	ALA B 185	-71.123	-36.856	-30.281	1.00	28.64	C
ATOM	10782	CB	ALA B 185	-71.689	-38.233	-29.986	1.00	28.44	C
ATOM	10786	C	ALA B 185	-70.017	-36.965	-31.330	1.00	28.08	C
ATOM	10787	O	ALA B 185	-68.839	-37.043	-30.992	1.00	28.26	O
ATOM	10789	N	LEU B 186	-70.394	-36.984	-32.602	1.00	27.40	N
ATOM	10790	CA	LEU B 186	-69.412	-37.101	-33.674	1.00	26.84	C
ATOM	10792	CB	LEU B 186	-70.088	-37.481	-35.000	1.00	26.75	C
ATOM	10795	CG	LEU B 186	-70.085	-38.983	-35.320	1.00	26.70	C
ATOM	10797	CD1	LEU B 186	-70.214	-39.868	-34.075	1.00	27.07	C
ATOM	10801	CD2	LEU B 186	-71.179	-39.306	-36.303	1.00	26.29	C
ATOM	10805	C	LEU B 186	-68.594	-35.822	-33.815	1.00	26.45	C
ATOM	10806	O	LEU B 186	-67.449	-35.875	-34.237	1.00	26.48	O
ATOM	10808	N	GLU B 187	-69.186	-34.685	-33.454	1.00	25.93	N
ATOM	10809	CA	GLU B 187	-68.479	-33.406	-33.403	1.00	25.57	C

ATOM	10811	CB	GLU B 187	-69.447	-32.310	-32.962	1.00	25.61	C
ATOM	10814	CG	GLU B 187	-69.035	-30.899	-33.325	1.00	26.29	C
ATOM	10817	CD	GLU B 187	-69.930	-29.857	-32.671	1.00	27.10	C
ATOM	10818	OE1	GLU B 187	-70.312	-30.057	-31.487	1.00	26.51	O
ATOM	10819	OE2	GLU B 187	-70.251	-28.846	-33.344	1.00	27.54	O
ATOM	10820	C	GLU B 187	-67.307	-33.490	-32.418	1.00	25.13	C
ATOM	10821	O	GLU B 187	-66.155	-33.165	-32.749	1.00	24.85	O
ATOM	10823	N	LEU B 188	-67.625	-33.938	-31.204	1.00	24.58	N
ATOM	10824	CA	LEU B 188	-66.644	-34.137	-30.148	1.00	24.06	C
ATOM	10826	CB	LEU B 188	-66.343	-32.817	-29.451	1.00	24.12	C
ATOM	10829	CG	LEU B 188	-65.042	-32.714	-28.670	1.00	23.61	C
ATOM	10831	CD1	LEU B 188	-63.895	-32.872	-29.629	1.00	23.12	C
ATOM	10835	CD2	LEU B 188	-64.976	-31.368	-27.964	1.00	23.18	C
ATOM	10839	C	LEU B 188	-67.248	-35.087	-29.140	1.00	23.70	C
ATOM	10840	O	LEU B 188	-68.392	-34.901	-28.743	1.00	23.68	O
ATOM	10842	N	PRO B 189	-66.493	-36.106	-28.713	1.00	23.34	N
ATOM	10843	CA	PRO B 189	-67.031	-36.983	-27.692	1.00	23.08	C
ATOM	10845	CB	PRO B 189	-66.018	-38.115	-27.638	1.00	22.96	C
ATOM	10848	CG	PRO B 189	-64.743	-37.452	-27.940	1.00	23.22	C
ATOM	10851	CD	PRO B 189	-65.059	-36.350	-28.929	1.00	23.54	C
ATOM	10854	C	PRO B 189	-67.053	-36.229	-26.387	1.00	22.86	C
ATOM	10855	O	PRO B 189	-66.284	-35.285	-26.215	1.00	23.03	O
ATOM	10856	N	LEU B 190	-67.912	-36.628	-25.465	1.00	22.56	N
ATOM	10857	CA	LEU B 190	-68.152	-35.773	-24.323	1.00	22.41	C
ATOM	10859	CB	LEU B 190	-69.567	-35.970	-23.766	1.00	23.04	C
ATOM	10862	CG	LEU B 190	-69.853	-37.136	-22.848	1.00	23.41	C
ATOM	10864	CD1	LEU B 190	-69.303	-36.733	-21.481	1.00	25.02	C
ATOM	10868	CD2	LEU B 190	-71.342	-37.411	-22.820	1.00	21.69	C
ATOM	10872	C	LEU B 190	-67.065	-35.899	-23.266	1.00	21.56	C
ATOM	10873	O	LEU B 190	-66.860	-34.976	-22.477	1.00	21.67	O
ATOM	10875	N	HIS B 191	-66.320	-36.998	-23.283	1.00	20.39	N
ATOM	10876	CA	HIS B 191	-65.089	-37.037	-22.490	1.00	19.40	C
ATOM	10878	CB	HIS B 191	-64.399	-38.393	-22.597	1.00	19.19	C
ATOM	10881	CG	HIS B 191	-63.222	-38.530	-21.689	1.00	19.32	C
ATOM	10882	ND1	HIS B 191	-63.347	-38.562	-20.317	1.00	20.95	N
ATOM	10884	CE1	HIS B 191	-62.145	-38.659	-19.771	1.00	21.15	C
ATOM	10886	NE2	HIS B 191	-61.245	-38.689	-20.741	1.00	19.36	N
ATOM	10888	CD2	HIS B 191	-61.893	-38.604	-21.948	1.00	19.49	C
ATOM	10890	C	HIS B 191	-64.094	-35.913	-22.841	1.00	18.62	C
ATOM	10891	O	HIS B 191	-63.150	-35.677	-22.085	1.00	18.59	O
ATOM	10893	N	ARG B 192	-64.307	-35.216	-23.957	1.00	17.72	N
ATOM	10894	CA	ARG B 192	-63.394	-34.150	-24.394	1.00	17.57	C
ATOM	10896	CB	ARG B 192	-62.838	-34.504	-25.773	1.00	17.81	C
ATOM	10899	CG	ARG B 192	-61.971	-35.736	-25.781	1.00	18.59	C
ATOM	10902	CD	ARG B 192	-61.484	-36.067	-27.181	1.00	20.10	C
ATOM	10905	NE	ARG B 192	-60.462	-37.116	-27.165	1.00	21.62	N
ATOM	10907	CZ	ARG B 192	-59.919	-37.648	-28.254	1.00	22.77	C
ATOM	10908	NH1	ARG B 192	-60.304	-37.248	-29.454	1.00	24.20	N
ATOM	10911	NH2	ARG B 192	-58.987	-38.581	-28.149	1.00	23.37	N
ATOM	10914	C	ARG B 192	-63.983	-32.726	-24.440	1.00	17.01	C
ATOM	10915	O	ARG B 192	-63.242	-31.739	-24.526	1.00	15.78	O
ATOM	10917	N	ARG B 193	-65.309	-32.634	-24.392	1.00	17.01	N
ATOM	10918	CA	ARG B 193	-66.015	-31.350	-24.453	1.00	16.99	C
ATOM	10920	CB	ARG B 193	-67.476	-31.582	-24.892	1.00	17.29	C
ATOM	10923	CG	ARG B 193	-68.192	-30.318	-25.377	1.00	18.42	C
ATOM	10926	CD	ARG B 193	-69.664	-30.559	-25.646	1.00	19.63	C
ATOM	10929	NE	ARG B 193	-69.912	-31.592	-26.648	1.00	21.16	N
ATOM	10931	CZ	ARG B 193	-69.868	-31.398	-27.971	1.00	23.56	C
ATOM	10932	NH1	ARG B 193	-69.575	-30.210	-28.496	1.00	23.37	N
ATOM	10935	NH2	ARG B 193	-70.112	-32.413	-28.789	1.00	25.37	N

ATOM	10938	C	ARG B 193	-65.974	-30.658	-23.087	1.00	16.46	C
ATOM	10939	O	ARG B 193	-66.040	-31.323	-22.034	1.00	16.47	O
ATOM	10941	N	THR B 194	-65.864	-29.335	-23.078	1.00	15.84	N
ATOM	10942	CA	THR B 194	-65.886	-28.629	-21.799	1.00	15.67	C
ATOM	10944	CB	THR B 194	-65.354	-27.199	-21.884	1.00	15.48	C
ATOM	10946	OG1	THR B 194	-66.077	-26.478	-22.882	1.00	15.55	O
ATOM	10948	CG2	THR B 194	-63.877	-27.189	-22.211	1.00	14.65	C
ATOM	10952	C	THR B 194	-67.310	-28.613	-21.270	1.00	15.79	C
ATOM	10953	O	THR B 194	-68.264	-28.652	-22.030	1.00	15.73	O
ATOM	10955	N	GLN B 195	-67.444	-28.557	-19.957	1.00	16.20	N
ATOM	10956	CA	GLN B 195	-68.737	-28.711	-19.328	1.00	16.84	C
ATOM	10958	CB	GLN B 195	-68.601	-28.781	-17.816	1.00	17.07	C
ATOM	10961	CG	GLN B 195	-69.921	-28.440	-17.153	1.00	18.86	C
ATOM	10964	CD	GLN B 195	-69.999	-28.875	-15.746	1.00	20.95	C
ATOM	10965	OE1	GLN B 195	-69.099	-29.511	-15.245	1.00	24.63	O
ATOM	10966	NE2	GLN B 195	-71.071	-28.529	-15.082	1.00	22.36	N
ATOM	10969	C	GLN B 195	-69.757	-27.623	-19.669	1.00	16.93	C
ATOM	10970	O	GLN B 195	-70.858	-27.931	-20.136	1.00	16.98	O
ATOM	10972	N	ARG B 196	-69.430	-26.366	-19.386	1.00	16.92	N
ATOM	10973	CA	ARG B 196	-70.358	-25.301	-19.685	1.00	17.06	C
ATOM	10975	CB	ARG B 196	-69.719	-23.932	-19.459	1.00	17.00	C
ATOM	10978	CG	ARG B 196	-70.095	-23.289	-18.120	1.00	17.09	C
ATOM	10981	CD	ARG B 196	-70.283	-24.320	-16.995	1.00	17.56	C
ATOM	10984	NE	ARG B 196	-71.452	-24.043	-16.146	1.00	17.64	N
ATOM	10986	CZ	ARG B 196	-72.289	-24.965	-15.664	1.00	17.93	C
ATOM	10987	NH1	ARG B 196	-72.136	-26.244	-15.966	1.00	17.72	N
ATOM	10990	NH2	ARG B 196	-73.305	-24.611	-14.878	1.00	18.27	N
ATOM	10993	C	ARG B 196	-70.881	-25.481	-21.098	1.00	17.42	C
ATOM	10994	O	ARG B 196	-72.079	-25.519	-21.317	1.00	17.66	O
ATOM	10996	N	LEU B 197	-69.984	-25.676	-22.044	1.00	18.04	N
ATOM	10997	CA	LEU B 197	-70.379	-25.925	-23.425	1.00	18.51	C
ATOM	10999	CB	LEU B 197	-69.133	-26.085	-24.289	1.00	18.52	C
ATOM	11002	CG	LEU B 197	-68.998	-25.185	-25.508	1.00	18.06	C
ATOM	11004	CD1	LEU B 197	-68.209	-23.932	-25.205	1.00	15.52	C
ATOM	11008	CD2	LEU B 197	-68.293	-26.002	-26.574	1.00	20.21	C
ATOM	11012	C	LEU B 197	-71.290	-27.163	-23.580	1.00	19.14	C
ATOM	11013	O	LEU B 197	-72.226	-27.150	-24.372	1.00	19.52	O
ATOM	11015	N	GLU B 198	-71.030	-28.235	-22.842	1.00	19.72	N
ATOM	11016	CA	GLU B 198	-71.918	-29.396	-22.909	1.00	20.43	C
ATOM	11018	CB	GLU B 198	-71.272	-30.631	-22.270	1.00	20.65	C
ATOM	11021	CG	GLU B 198	-72.176	-31.880	-22.172	1.00	22.51	C
ATOM	11024	CD	GLU B 198	-72.635	-32.460	-23.522	1.00	25.57	C
ATOM	11025	OE1	GLU B 198	-72.315	-31.896	-24.581	1.00	28.90	O
ATOM	11026	OE2	GLU B 198	-73.326	-33.506	-23.540	1.00	28.09	O
ATOM	11027	C	GLU B 198	-73.271	-29.101	-22.259	1.00	20.72	C
ATOM	11028	O	GLU B 198	-74.281	-29.693	-22.644	1.00	20.82	O
ATOM	11030	N	ALA B 199	-73.282	-28.199	-21.270	1.00	20.86	N
ATOM	11031	CA	ALA B 199	-74.508	-27.814	-20.565	1.00	20.66	C
ATOM	11033	CB	ALA B 199	-74.176	-27.013	-19.326	1.00	20.36	C
ATOM	11037	C	ALA B 199	-75.472	-27.033	-21.461	1.00	20.84	C
ATOM	11038	O	ALA B 199	-76.644	-27.404	-21.590	1.00	20.84	O
ATOM	11040	N	VAL B 200	-75.002	-25.961	-22.092	1.00	21.12	N
ATOM	11041	CA	VAL B 200	-75.903	-25.177	-22.942	1.00	21.60	C
ATOM	11043	CB	VAL B 200	-75.215	-24.042	-23.701	1.00	21.30	C
ATOM	11045	CG1	VAL B 200	-74.660	-23.055	-22.738	1.00	21.26	C
ATOM	11049	CG2	VAL B 200	-74.135	-24.575	-24.626	1.00	21.28	C
ATOM	11053	C	VAL B 200	-76.603	-26.077	-23.944	1.00	22.32	C
ATOM	11054	O	VAL B 200	-77.791	-25.908	-24.218	1.00	22.32	O
ATOM	11056	N	TRP B 201	-75.872	-27.047	-24.472	1.00	23.16	N
ATOM	11057	CA	TRP B 201	-76.441	-27.946	-25.454	1.00	24.02	C

ATOM	11059	CB	TRP	B	201	-75.348	-28.718	-26.195	1.00	24.36	C
ATOM	11062	CG	TRP	B	201	-75.898	-29.510	-27.331	1.00	24.84	C
ATOM	11063	CD1	TRP	B	201	-76.060	-29.092	-28.613	1.00	25.41	C
ATOM	11065	NE1	TRP	B	201	-76.599	-30.094	-29.374	1.00	25.61	N
ATOM	11067	CE2	TRP	B	201	-76.807	-31.184	-28.579	1.00	25.34	C
ATOM	11068	CD2	TRP	B	201	-76.374	-30.847	-27.279	1.00	25.56	C
ATOM	11069	CE3	TRP	B	201	-76.466	-31.800	-26.264	1.00	27.13	C
ATOM	11071	CZ3	TRP	B	201	-76.989	-33.046	-26.574	1.00	28.66	C
ATOM	11073	CH2	TRP	B	201	-77.419	-33.347	-27.885	1.00	28.01	C
ATOM	11075	CZ2	TRP	B	201	-77.332	-32.425	-28.894	1.00	26.09	C
ATOM	11077	C	TRP	B	201	-77.411	-28.938	-24.833	1.00	24.22	C
ATOM	11078	O	TRP	B	201	-78.473	-29.172	-25.383	1.00	24.83	O
ATOM	11080	N	SER	B	202	-77.034	-29.541	-23.715	1.00	24.39	N
ATOM	11081	CA	SER	B	202	-77.818	-30.626	-23.144	1.00	24.50	C
ATOM	11083	CB	SER	B	202	-76.968	-31.455	-22.184	1.00	24.42	C
ATOM	11086	OG	SER	B	202	-75.848	-32.021	-22.848	1.00	23.80	O
ATOM	11088	C	SER	B	202	-79.067	-30.099	-22.441	1.00	24.98	C
ATOM	11089	O	SER	B	202	-80.057	-30.816	-22.324	1.00	24.73	O
ATOM	11091	N	ILE	B	203	-79.032	-28.849	-21.982	1.00	25.60	N
ATOM	11092	CA	ILE	B	203	-80.219	-28.255	-21.372	1.00	25.98	C
ATOM	11094	CB	ILE	B	203	-79.908	-26.967	-20.591	1.00	25.85	C
ATOM	11096	CG1	ILE	B	203	-79.086	-27.326	-19.348	1.00	25.76	C
ATOM	11099	CD1	ILE	B	203	-78.779	-26.159	-18.431	1.00	26.28	C
ATOM	11103	CG2	ILE	B	203	-81.212	-26.253	-20.206	1.00	24.51	C
ATOM	11107	C	ILE	B	203	-81.280	-28.016	-22.443	1.00	26.67	C
ATOM	11108	O	ILE	B	203	-82.428	-28.398	-22.260	1.00	26.96	O
ATOM	11110	N	GLU	B	204	-80.879	-27.412	-23.563	1.00	27.30	N
ATOM	11111	CA	GLU	B	204	-81.758	-27.205	-24.720	1.00	27.49	C
ATOM	11113	CB	GLU	B	204	-81.036	-26.384	-25.804	1.00	27.67	C
ATOM	11116	CG	GLU	B	204	-81.849	-26.077	-27.073	1.00	28.50	C
ATOM	11119	CD	GLU	B	204	-82.987	-25.072	-26.855	1.00	29.81	C
ATOM	11120	OE1	GLU	B	204	-83.158	-24.600	-25.711	1.00	31.40	O
ATOM	11121	OE2	GLU	B	204	-83.710	-24.750	-27.832	1.00	29.57	O
ATOM	11122	C	GLU	B	204	-82.201	-28.542	-25.287	1.00	27.52	C
ATOM	11123	O	GLU	B	204	-83.321	-28.667	-25.738	1.00	27.67	O
ATOM	11125	N	ALA	B	205	-81.326	-29.539	-25.266	1.00	27.74	N
ATOM	11126	CA	ALA	B	205	-81.684	-30.872	-25.728	1.00	28.09	C
ATOM	11128	CB	ALA	B	205	-80.472	-31.789	-25.726	1.00	27.75	C
ATOM	11132	C	ALA	B	205	-82.771	-31.448	-24.839	1.00	28.81	C
ATOM	11133	O	ALA	B	205	-83.811	-31.877	-25.325	1.00	28.90	O
ATOM	11135	N	TYR	B	206	-82.517	-31.433	-23.531	1.00	29.72	N
ATOM	11136	CA	TYR	B	206	-83.408	-32.020	-22.527	1.00	30.30	C
ATOM	11138	CB	TYR	B	206	-82.760	-31.915	-21.149	1.00	30.24	C
ATOM	11141	CG	TYR	B	206	-83.276	-32.889	-20.127	1.00	29.79	C
ATOM	11142	CD1	TYR	B	206	-83.055	-34.245	-20.286	1.00	30.71	C
ATOM	11144	CE1	TYR	B	206	-83.498	-35.161	-19.360	1.00	30.78	C
ATOM	11146	CZ	TYR	B	206	-84.159	-34.731	-18.238	1.00	30.03	C
ATOM	11147	OH	TYR	B	206	-84.582	-35.684	-17.343	1.00	30.51	O
ATOM	11149	CE2	TYR	B	206	-84.394	-33.379	-18.042	1.00	29.42	C
ATOM	11151	CD2	TYR	B	206	-83.945	-32.463	-18.989	1.00	29.05	C
ATOM	11153	C	TYR	B	206	-84.744	-31.308	-22.467	1.00	31.19	C
ATOM	11154	O	TYR	B	206	-85.790	-31.930	-22.311	1.00	31.44	O
ATOM	11156	N	ARG	B	207	-84.683	-29.988	-22.558	1.00	32.21	N
ATOM	11157	CA	ARG	B	207	-85.853	-29.128	-22.521	1.00	33.04	C
ATOM	11159	CB	ARG	B	207	-85.394	-27.700	-22.843	1.00	32.86	C
ATOM	11162	CG	ARG	B	207	-86.441	-26.637	-22.846	1.00	32.82	C
ATOM	11165	CD	ARG	B	207	-85.905	-25.360	-23.462	1.00	32.76	C
ATOM	11168	NE	ARG	B	207	-85.148	-24.561	-22.504	1.00	32.80	N
ATOM	11170	CZ	ARG	B	207	-85.676	-23.709	-21.628	1.00	33.40	C
ATOM	11171	NH1	ARG	B	207	-86.989	-23.526	-21.555	1.00	34.29	N

ATOM	11174	NH2 ARG B 207	-84.884	-23.030	-20.808	1.00	33.89	N
ATOM	11177	C ARG B 207	-86.971	-29.606	-23.478	1.00	34.14	C
ATOM	11178	O ARG B 207	-88.144	-29.521	-23.128	1.00	34.34	O
ATOM	11180	N LYS B 208	-86.607	-30.137	-24.652	1.00	35.32	N
ATOM	11181	CA LYS B 208	-87.581	-30.535	-25.688	1.00	36.16	C
ATOM	11183	CB LYS B 208	-86.960	-30.468	-27.090	1.00	36.01	C
ATOM	11186	CG LYS B 208	-86.126	-29.234	-27.368	1.00	35.85	C
ATOM	11189	CD LYS B 208	-85.934	-29.023	-28.861	1.00	35.98	C
ATOM	11192	CE LYS B 208	-84.774	-28.086	-29.186	1.00	35.63	C
ATOM	11195	NZ LYS B 208	-83.566	-28.843	-29.596	1.00	35.17	N
ATOM	11199	C LYS B 208	-88.166	-31.937	-25.490	1.00	37.25	C
ATOM	11200	O LYS B 208	-89.212	-32.250	-26.054	1.00	37.50	O
ATOM	11202	N LYS B 209	-87.484	-32.788	-24.730	1.00	38.58	N
ATOM	11203	CA LYS B 209	-88.026	-34.106	-24.383	1.00	39.89	C
ATOM	11205	CB LYS B 209	-87.013	-34.923	-23.578	1.00	40.23	C
ATOM	11208	CG LYS B 209	-85.914	-35.604	-24.373	1.00	41.47	C
ATOM	11211	CD LYS B 209	-84.858	-36.178	-23.417	1.00	43.71	C
ATOM	11214	CE LYS B 209	-84.448	-37.612	-23.764	1.00	45.15	C
ATOM	11217	NZ LYS B 209	-85.390	-38.621	-23.163	1.00	45.73	N
ATOM	11221	C LYS B 209	-89.288	-33.950	-23.530	1.00	40.51	C
ATOM	11222	O LYS B 209	-89.230	-33.367	-22.441	1.00	40.84	O
ATOM	11224	N GLU B 210	-90.418	-34.477	-24.000	1.00	40.92	N
ATOM	11225	CA GLU B 210	-91.652	-34.386	-23.223	1.00	41.27	C
ATOM	11227	CB GLU B 210	-92.855	-34.935	-23.998	1.00	41.77	C
ATOM	11230	CG GLU B 210	-92.822	-36.446	-24.265	1.00	43.47	C
ATOM	11233	CD GLU B 210	-94.049	-36.922	-25.032	1.00	45.23	C
ATOM	11234	OE1 GLU B 210	-95.188	-36.590	-24.612	1.00	45.02	O
ATOM	11235	OE2 GLU B 210	-93.863	-37.629	-26.053	1.00	46.74	O
ATOM	11236	C GLU B 210	-91.494	-35.121	-21.895	1.00	40.71	C
ATOM	11237	O GLU B 210	-91.996	-34.667	-20.864	1.00	40.99	O
ATOM	11239	N ASP B 211	-90.773	-36.240	-21.924	1.00	39.76	N
ATOM	11240	CA ASP B 211	-90.533	-37.050	-20.721	1.00	38.96	C
ATOM	11242	CB ASP B 211	-90.151	-38.476	-21.123	1.00	39.10	C
ATOM	11245	CG ASP B 211	-89.101	-38.499	-22.212	1.00	40.26	C
ATOM	11246	OD1 ASP B 211	-89.327	-37.827	-23.254	1.00	41.25	O
ATOM	11247	OD2 ASP B 211	-88.055	-39.160	-22.020	1.00	41.66	O
ATOM	11248	C ASP B 211	-89.440	-36.468	-19.826	1.00	37.75	C
ATOM	11249	O ASP B 211	-89.009	-37.131	-18.892	1.00	37.62	O
ATOM	11251	N ALA B 212	-88.989	-35.245	-20.115	1.00	36.47	N
ATOM	11252	CA ALA B 212	-87.941	-34.587	-19.335	1.00	35.36	C
ATOM	11254	CB ALA B 212	-87.516	-33.294	-20.011	1.00	35.23	C
ATOM	11258	C ALA B 212	-88.423	-34.307	-17.920	1.00	34.40	C
ATOM	11259	O ALA B 212	-89.559	-33.870	-17.728	1.00	34.42	O
ATOM	11261	N ASN B 213	-87.565	-34.574	-16.937	1.00	33.24	N
ATOM	11262	CA ASN B 213	-87.890	-34.333	-15.533	1.00	32.67	C
ATOM	11264	CB ASN B 213	-86.840	-34.986	-14.623	1.00	32.58	C
ATOM	11267	CG ASN B 213	-87.204	-34.912	-13.152	1.00	32.67	C
ATOM	11268	OD1 ASN B 213	-87.891	-33.994	-12.722	1.00	32.62	O
ATOM	11269	ND2 ASN B 213	-86.734	-35.882	-12.371	1.00	33.18	N
ATOM	11272	C ASN B 213	-87.990	-32.825	-15.266	1.00	32.19	C
ATOM	11273	O ASN B 213	-87.010	-32.092	-15.376	1.00	32.38	O
ATOM	11275	N GLN B 214	-89.182	-32.356	-14.923	1.00	31.48	N
ATOM	11276	CA GLN B 214	-89.389	-30.924	-14.756	1.00	30.85	C
ATOM	11278	CB GLN B 214	-90.889	-30.581	-14.693	1.00	30.82	C
ATOM	11281	CG GLN B 214	-91.684	-30.869	-16.002	1.00	30.80	C
ATOM	11284	CD GLN B 214	-90.990	-30.378	-17.291	1.00	29.69	C
ATOM	11285	OE1 GLN B 214	-90.893	-29.174	-17.543	1.00	29.07	O
ATOM	11286	NE2 GLN B 214	-90.526	-31.321	-18.113	1.00	27.71	N
ATOM	11289	C GLN B 214	-88.638	-30.369	-13.543	1.00	30.37	C
ATOM	11290	O GLN B 214	-88.263	-29.200	-13.533	1.00	30.80	O

ATOM	11292	N	VAL B 215	-88.393	-31.200	-12.532	1.00	29.54	N
ATOM	11293	CA	VAL B 215	-87.566	-30.779	-11.395	1.00	28.56	C
ATOM	11295	CB	VAL B 215	-87.564	-31.818	-10.260	1.00	28.64	C
ATOM	11297	CG1	VAL B 215	-86.758	-31.299	-9.075	1.00	27.94	C
ATOM	11301	CG2	VAL B 215	-88.999	-32.167	-9.856	1.00	28.26	C
ATOM	11305	C	VAL B 215	-86.124	-30.545	-11.845	1.00	27.63	C
ATOM	11306	O	VAL B 215	-85.566	-29.472	-11.615	1.00	27.58	O
ATOM	11308	N	LEU B 216	-85.546	-31.548	-12.504	1.00	26.37	N
ATOM	11309	CA	LEU B 216	-84.143	-31.509	-12.930	1.00	25.47	C
ATOM	11311	CB	LEU B 216	-83.714	-32.866	-13.490	1.00	25.46	C
ATOM	11314	CG	LEU B 216	-82.274	-32.995	-13.979	1.00	24.97	C
ATOM	11316	CD1	LEU B 216	-81.289	-32.640	-12.883	1.00	24.62	C
ATOM	11320	CD2	LEU B 216	-82.042	-34.407	-14.472	1.00	24.55	C
ATOM	11324	C	LEU B 216	-83.861	-30.447	-13.972	1.00	24.68	C
ATOM	11325	O	LEU B 216	-82.819	-29.820	-13.929	1.00	24.91	O
ATOM	11327	N	LEU B 217	-84.776	-30.269	-14.916	1.00	23.79	N
ATOM	11328	CA	LEU B 217	-84.625	-29.260	-15.960	1.00	23.13	C
ATOM	11330	CB	LEU B 217	-85.765	-29.384	-16.972	1.00	23.36	C
ATOM	11333	CG	LEU B 217	-85.808	-28.374	-18.123	1.00	23.18	C
ATOM	11335	CD1	LEU B 217	-84.462	-28.268	-18.796	1.00	22.69	C
ATOM	11339	CD2	LEU B 217	-86.861	-28.786	-19.126	1.00	23.38	C
ATOM	11343	C	LEU B 217	-84.631	-27.853	-15.383	1.00	22.51	C
ATOM	11344	O	LEU B 217	-83.903	-26.973	-15.855	1.00	22.08	O
ATOM	11346	N	GLU B 218	-85.481	-27.647	-14.375	1.00	21.84	N
ATOM	11347	CA	GLU B 218	-85.619	-26.339	-13.725	1.00	21.26	C
ATOM	11349	CB	GLU B 218	-86.813	-26.341	-12.768	1.00	21.27	C
ATOM	11352	CG	GLU B 218	-87.206	-24.956	-12.246	1.00	21.35	C
ATOM	11355	CD	GLU B 218	-88.501	-24.964	-11.441	1.00	20.89	C
ATOM	11356	OE1	GLU B 218	-89.097	-26.052	-11.263	1.00	20.47	O
ATOM	11357	OE2	GLU B 218	-88.915	-23.874	-10.986	1.00	20.15	O
ATOM	11358	C	GLU B 218	-84.348	-26.001	-12.969	1.00	20.58	C
ATOM	11359	O	GLU B 218	-83.877	-24.865	-12.982	1.00	20.31	O
ATOM	11361	N	LEU B 219	-83.802	-27.015	-12.312	1.00	19.87	N
ATOM	11362	CA	LEU B 219	-82.559	-26.885	-11.592	1.00	19.12	C
ATOM	11364	CB	LEU B 219	-82.310	-28.132	-10.746	1.00	18.58	C
ATOM	11367	CG	LEU B 219	-81.159	-28.054	-9.762	1.00	17.07	C
ATOM	11369	CD1	LEU B 219	-81.321	-26.847	-8.872	1.00	16.10	C
ATOM	11373	CD2	LEU B 219	-81.107	-29.311	-8.954	1.00	15.22	C
ATOM	11377	C	LEU B 219	-81.443	-26.675	-12.599	1.00	19.16	C
ATOM	11378	O	LEU B 219	-80.684	-25.721	-12.478	1.00	19.42	O
ATOM	11380	N	ALA B 220	-81.368	-27.546	-13.606	1.00	19.04	N
ATOM	11381	CA	ALA B 220	-80.356	-27.431	-14.669	1.00	19.04	C
ATOM	11383	CB	ALA B 220	-80.648	-28.398	-15.817	1.00	18.64	C
ATOM	11387	C	ALA B 220	-80.257	-26.004	-15.196	1.00	19.06	C
ATOM	11388	O	ALA B 220	-79.159	-25.472	-15.321	1.00	19.08	O
ATOM	11390	N	ILE B 221	-81.410	-25.393	-15.475	1.00	19.19	N
ATOM	11391	CA	ILE B 221	-81.473	-24.016	-15.963	1.00	19.22	C
ATOM	11393	CB	ILE B 221	-82.898	-23.638	-16.423	1.00	19.03	C
ATOM	11395	CG1	ILE B 221	-83.229	-24.323	-17.755	1.00	18.52	C
ATOM	11398	CD1	ILE B 221	-84.687	-24.645	-17.928	1.00	17.52	C
ATOM	11402	CG2	ILE B 221	-83.024	-22.121	-16.564	1.00	18.03	C
ATOM	11406	C	ILE B 221	-81.030	-23.013	-14.906	1.00	19.62	C
ATOM	11407	O	ILE B 221	-80.137	-22.190	-15.136	1.00	19.53	O
ATOM	11409	N	LEU B 222	-81.666	-23.096	-13.745	1.00	20.27	N
ATOM	11410	CA	LEU B 222	-81.412	-22.161	-12.654	1.00	20.63	C
ATOM	11412	CB	LEU B 222	-82.203	-22.561	-11.412	1.00	20.50	C
ATOM	11415	CG	LEU B 222	-82.195	-21.523	-10.296	1.00	20.76	C
ATOM	11417	CD1	LEU B 222	-83.439	-21.675	-9.431	1.00	21.02	C
ATOM	11421	CD2	LEU B 222	-80.924	-21.625	-9.453	1.00	20.21	C
ATOM	11425	C	LEU B 222	-79.925	-22.089	-12.341	1.00	20.98	C

ATOM	11426	O	LEU B 222	-79.337	-21.011	-12.381	1.00	21.40	O
ATOM	11428	N	ASP B 223	-79.319	-23.240	-12.064	1.00	21.24	N
ATOM	11429	CA	ASP B 223	-77.907	-23.298	-11.709	1.00	21.48	C
ATOM	11431	CB	ASP B 223	-77.509	-24.724	-11.314	1.00	21.66	C
ATOM	11434	CG	ASP B 223	-76.168	-24.792	-10.589	1.00	22.70	C
ATOM	11435	OD1	ASP B 223	-75.095	-24.775	-11.244	1.00	24.25	O
ATOM	11436	OD2	ASP B 223	-76.190	-24.893	-9.349	1.00	24.85	O
ATOM	11437	C	ASP B 223	-77.017	-22.768	-12.835	1.00	21.64	C
ATOM	11438	O	ASP B 223	-76.035	-22.104	-12.548	1.00	21.66	O
ATOM	11440	N	TYR B 224	-77.361	-23.023	-14.100	1.00	22.07	N
ATOM	11441	CA	TYR B 224	-76.517	-22.561	-15.218	1.00	22.47	C
ATOM	11443	CB	TYR B 224	-76.980	-23.111	-16.589	1.00	22.35	C
ATOM	11446	CG	TYR B 224	-76.032	-22.733	-17.724	1.00	22.25	C
ATOM	11447	CD1	TYR B 224	-75.002	-23.576	-18.116	1.00	21.46	C
ATOM	11449	CE1	TYR B 224	-74.126	-23.216	-19.135	1.00	21.43	C
ATOM	11451	CZ	TYR B 224	-74.260	-21.991	-19.761	1.00	21.62	C
ATOM	11452	OH	TYR B 224	-73.387	-21.618	-20.760	1.00	20.46	O
ATOM	11454	CE2	TYR B 224	-75.270	-21.132	-19.387	1.00	22.04	C
ATOM	11456	CD2	TYR B 224	-76.144	-21.500	-18.373	1.00	22.73	C
ATOM	11458	C	TYR B 224	-76.414	-21.023	-15.262	1.00	23.06	C
ATOM	11459	O	TYR B 224	-75.323	-20.467	-15.507	1.00	22.89	O
ATOM	11461	N	ASN B 225	-77.543	-20.350	-15.027	1.00	23.66	N
ATOM	11462	CA	ASN B 225	-77.593	-18.884	-15.072	1.00	24.16	C
ATOM	11464	CB	ASN B 225	-79.040	-18.367	-15.107	1.00	24.12	C
ATOM	11467	CG	ASN B 225	-79.778	-18.761	-16.365	1.00	23.89	C
ATOM	11468	OD1	ASN B 225	-79.182	-18.877	-17.437	1.00	24.33	O
ATOM	11469	ND2	ASN B 225	-81.088	-18.965	-16.244	1.00	22.05	N
ATOM	11472	C	ASN B 225	-76.880	-18.265	-13.881	1.00	24.73	C
ATOM	11473	O	ASN B 225	-76.195	-17.247	-14.027	1.00	25.09	O
ATOM	11475	N	MET B 226	-77.067	-18.864	-12.703	1.00	25.15	N
ATOM	11476	CA	MET B 226	-76.423	-18.390	-11.479	1.00	25.53	C
ATOM	11478	CB	MET B 226	-76.806	-19.278	-10.282	1.00	26.11	C
ATOM	11481	CG	MET B 226	-75.905	-19.169	-9.023	1.00	27.59	C
ATOM	11484	SD	MET B 226	-75.676	-20.788	-8.214	1.00	30.67	S
ATOM	11485	CE	MET B 226	-77.308	-21.005	-7.476	1.00	29.78	C
ATOM	11489	C	MET B 226	-74.920	-18.380	-11.692	1.00	25.16	C
ATOM	11490	O	MET B 226	-74.258	-17.387	-11.391	1.00	25.10	O
ATOM	11492	N	ILE B 227	-74.388	-19.473	-12.239	1.00	24.88	N
ATOM	11493	CA	ILE B 227	-72.944	-19.564	-12.501	1.00	24.70	C
ATOM	11495	CB	ILE B 227	-72.476	-20.994	-12.882	1.00	24.45	C
ATOM	11497	CG1	ILE B 227	-72.656	-21.947	-11.695	1.00	23.80	C
ATOM	11500	CD1	ILE B 227	-72.094	-23.332	-11.885	1.00	21.97	C
ATOM	11504	CG2	ILE B 227	-71.030	-20.969	-13.266	1.00	24.90	C
ATOM	11508	C	ILE B 227	-72.529	-18.550	-13.566	1.00	24.45	C
ATOM	11509	O	ILE B 227	-71.578	-17.810	-13.373	1.00	24.42	O
ATOM	11511	N	GLN B 228	-73.264	-18.484	-14.668	1.00	24.31	N
ATOM	11512	CA	GLN B 228	-73.038	-17.414	-15.639	1.00	24.19	C
ATOM	11514	CB	GLN B 228	-74.143	-17.363	-16.688	1.00	24.02	C
ATOM	11517	CG	GLN B 228	-73.788	-16.488	-17.864	1.00	23.42	C
ATOM	11520	CD	GLN B 228	-74.807	-16.567	-18.983	1.00	23.70	C
ATOM	11521	OE1	GLN B 228	-75.790	-15.813	-19.015	1.00	24.07	O
ATOM	11522	NE2	GLN B 228	-74.564	-17.465	-19.927	1.00	23.49	N
ATOM	11525	C	GLN B 228	-72.913	-16.037	-14.987	1.00	24.37	C
ATOM	11526	O	GLN B 228	-72.145	-15.210	-15.470	1.00	24.45	O
ATOM	11528	N	SER B 229	-73.653	-15.776	-13.909	1.00	24.45	N
ATOM	11529	CA	SER B 229	-73.578	-14.459	-13.264	1.00	24.76	C
ATOM	11531	CB	SER B 229	-74.815	-14.163	-12.397	1.00	24.84	C
ATOM	11534	OG	SER B 229	-74.711	-14.718	-11.096	1.00	25.81	O
ATOM	11536	C	SER B 229	-72.270	-14.279	-12.469	1.00	24.64	C
ATOM	11537	O	SER B 229	-71.726	-13.172	-12.408	1.00	24.95	O

ATOM	11539	N	VAL B 230	-71.751	-15.350	-11.878	1.00	24.20	N
ATOM	11540	CA	VAL B 230	-70.415	-15.277	-11.299	1.00	23.99	C
ATOM	11542	CB	VAL B 230	-70.006	-16.581	-10.574	1.00	23.92	C
ATOM	11544	CG1	VAL B 230	-68.546	-16.510	-10.105	1.00	23.71	C
ATOM	11548	CG2	VAL B 230	-70.923	-16.853	-9.409	1.00	23.52	C
ATOM	11552	C	VAL B 230	-69.386	-14.966	-12.400	1.00	24.12	C
ATOM	11553	O	VAL B 230	-68.397	-14.291	-12.136	1.00	24.01	O
ATOM	11555	N	TYR B 231	-69.612	-15.461	-13.621	1.00	24.30	N
ATOM	11556	CA	TYR B 231	-68.672	-15.235	-14.728	1.00	24.53	C
ATOM	11558	CB	TYR B 231	-68.997	-16.104	-15.946	1.00	24.08	C
ATOM	11561	CG	TYR B 231	-68.892	-17.605	-15.754	1.00	23.29	C
ATOM	11562	CD1	TYR B 231	-68.173	-18.168	-14.706	1.00	23.08	C
ATOM	11564	CE1	TYR B 231	-68.073	-19.547	-14.565	1.00	22.01	C
ATOM	11566	CZ	TYR B 231	-68.685	-20.369	-15.479	1.00	21.16	C
ATOM	11567	OH	TYR B 231	-68.607	-21.746	-15.353	1.00	20.18	O
ATOM	11569	CE2	TYR B 231	-69.387	-19.825	-16.523	1.00	21.64	C
ATOM	11571	CD2	TYR B 231	-69.476	-18.464	-16.662	1.00	22.15	C
ATOM	11573	C	TYR B 231	-68.658	-13.772	-15.166	1.00	25.37	C
ATOM	11574	O	TYR B 231	-67.602	-13.221	-15.528	1.00	25.25	O
ATOM	11576	N	GLN B 232	-69.832	-13.148	-15.129	1.00	26.11	N
ATOM	11577	CA	GLN B 232	-69.978	-11.797	-15.622	1.00	26.68	C
ATOM	11579	CB	GLN B 232	-71.430	-11.541	-15.995	1.00	26.45	C
ATOM	11582	CG	GLN B 232	-71.883	-12.385	-17.176	1.00	25.67	C
ATOM	11585	CD	GLN B 232	-73.388	-12.394	-17.381	1.00	24.91	C
ATOM	11586	OE1	GLN B 232	-74.140	-11.761	-16.637	1.00	24.81	O
ATOM	11587	NE2	GLN B 232	-73.834	-13.113	-18.407	1.00	23.39	N
ATOM	11590	C	GLN B 232	-69.445	-10.802	-14.593	1.00	27.90	C
ATOM	11591	O	GLN B 232	-68.909	-9.752	-14.958	1.00	27.96	O
ATOM	11593	N	ARG B 233	-69.572	-11.145	-13.312	1.00	29.35	N
ATOM	11594	CA	ARG B 233	-68.908	-10.396	-12.241	1.00	30.59	C
ATOM	11596	CB	ARG B 233	-69.490	-10.775	-10.867	1.00	30.99	C
ATOM	11599	CG	ARG B 233	-68.824	-10.114	-9.661	1.00	32.72	C
ATOM	11602	CD	ARG B 233	-69.695	-10.210	-8.391	1.00	35.70	C
ATOM	11605	NE	ARG B 233	-70.162	-11.578	-8.086	1.00	38.29	N
ATOM	11607	CZ	ARG B 233	-71.403	-12.053	-8.276	1.00	39.65	C
ATOM	11608	NH1	ARG B 233	-72.367	-11.287	-8.787	1.00	40.23	N
ATOM	11611	NH2	ARG B 233	-71.687	-13.320	-7.949	1.00	39.61	N
ATOM	11614	C	ARG B 233	-67.390	-10.632	-12.303	1.00	31.10	C
ATOM	11615	O	ARG B 233	-66.615	-9.693	-12.181	1.00	31.19	O
ATOM	11617	N	ASP B 234	-66.954	-11.868	-12.517	1.00	31.91	N
ATOM	11618	CA	ASP B 234	-65.530	-12.102	-12.750	1.00	32.63	C
ATOM	11620	CB	ASP B 234	-65.220	-13.566	-13.090	1.00	32.53	C
ATOM	11623	CG	ASP B 234	-65.336	-14.493	-11.894	1.00	32.90	C
ATOM	11624	OD1	ASP B 234	-65.477	-14.028	-10.746	1.00	33.06	O
ATOM	11625	OD2	ASP B 234	-65.293	-15.713	-12.107	1.00	34.22	O
ATOM	11626	C	ASP B 234	-65.058	-11.206	-13.886	1.00	33.10	C
ATOM	11627	O	ASP B 234	-64.097	-10.458	-13.723	1.00	33.55	O
ATOM	11629	N	LEU B 235	-65.753	-11.256	-15.017	1.00	33.49	N
ATOM	11630	CA	LEU B 235	-65.311	-10.550	-16.216	1.00	34.06	C
ATOM	11632	CB	LEU B 235	-66.139	-10.996	-17.422	1.00	33.80	C
ATOM	11635	CG	LEU B 235	-65.769	-10.403	-18.775	1.00	32.13	C
ATOM	11637	CD1	LEU B 235	-64.303	-10.679	-19.050	1.00	30.94	C
ATOM	11641	CD2	LEU B 235	-66.670	-10.959	-19.870	1.00	29.93	C
ATOM	11645	C	LEU B 235	-65.345	-9.023	-16.103	1.00	35.40	C
ATOM	11646	O	LEU B 235	-64.500	-8.352	-16.682	1.00	35.61	O
ATOM	11648	N	ARG B 236	-66.327	-8.473	-15.386	1.00	36.92	N
ATOM	11649	CA	ARG B 236	-66.380	-7.019	-15.137	1.00	38.06	C
ATOM	11651	CB	ARG B 236	-67.650	-6.607	-14.367	1.00	38.43	C
ATOM	11654	CG	ARG B 236	-68.882	-6.321	-15.230	1.00	39.58	C
ATOM	11657	CD	ARG B 236	-69.958	-5.520	-14.459	1.00	40.93	C

ATOM	11660	NE	ARG B 236	-70.383	-6.149	-13.199	1.00	41.74	N
ATOM	11662	CZ	ARG B 236	-71.208	-7.195	-13.095	1.00	42.20	C
ATOM	11663	NH1	ARG B 236	-71.719	-7.784	-14.175	1.00	41.61	N
ATOM	11666	NH2	ARG B 236	-71.517	-7.672	-11.892	1.00	42.99	N
ATOM	11669	C	ARG B 236	-65.170	-6.550	-14.343	1.00	38.57	C
ATOM	11670	O	ARG B 236	-64.593	-5.524	-14.665	1.00	38.61	O
ATOM	11672	N	GLU B 237	-64.822	-7.292	-13.291	1.00	39.42	N
ATOM	11673	CA	GLU B 237	-63.625	-7.017	-12.480	1.00	40.18	C
ATOM	11675	CB	GLU B 237	-63.471	-8.033	-11.317	1.00	40.67	C
ATOM	11678	CG	GLU B 237	-64.336	-7.750	-10.057	1.00	42.61	C
ATOM	11681	CD	GLU B 237	-64.258	-8.858	-8.961	1.00	45.07	C
ATOM	11682	OE1	GLU B 237	-64.051	-10.066	-9.282	1.00	46.74	O
ATOM	11683	OE2	GLU B 237	-64.436	-8.511	-7.764	1.00	45.87	O
ATOM	11684	C	GLU B 237	-62.352	-7.024	-13.344	1.00	39.93	C
ATOM	11685	O	GLU B 237	-61.593	-6.050	-13.332	1.00	39.96	O
ATOM	11687	N	THR B 238	-62.123	-8.104	-14.095	1.00	39.56	N
ATOM	11688	CA	THR B 238	-60.900	-8.200	-14.895	1.00	39.42	C
ATOM	11690	CB	THR B 238	-60.492	-9.666	-15.268	1.00	39.39	C
ATOM	11692	OG1	THR B 238	-61.222	-10.123	-16.413	1.00	39.37	O
ATOM	11694	CG2	THR B 238	-60.683	-10.612	-14.095	1.00	38.73	C
ATOM	11698	C	THR B 238	-60.961	-7.326	-16.156	1.00	39.48	C
ATOM	11699	O	THR B 238	-59.930	-7.062	-16.761	1.00	39.48	O
ATOM	11701	N	SER B 239	-62.149	-6.871	-16.549	1.00	39.62	N
ATOM	11702	CA	SER B 239	-62.264	-5.882	-17.631	1.00	39.63	C
ATOM	11704	CB	SER B 239	-63.689	-5.810	-18.172	1.00	39.56	C
ATOM	11707	OG	SER B 239	-63.945	-6.917	-19.010	1.00	38.96	O
ATOM	11709	C	SER B 239	-61.796	-4.496	-17.181	1.00	39.99	C
ATOM	11710	O	SER B 239	-61.108	-3.807	-17.933	1.00	40.10	O
ATOM	11712	N	ARG B 240	-62.168	-4.088	-15.965	1.00	40.35	N
ATOM	11713	CA	ARG B 240	-61.624	-2.864	-15.363	1.00	40.68	C
ATOM	11715	CB	ARG B 240	-62.025	-2.712	-13.881	1.00	41.26	C
ATOM	11718	CG	ARG B 240	-63.231	-1.782	-13.621	1.00	43.86	C
ATOM	11721	CD	ARG B 240	-63.334	-1.375	-12.130	1.00	47.07	C
ATOM	11724	NE	ARG B 240	-63.365	-2.537	-11.222	1.00	50.14	N
ATOM	11726	CZ	ARG B 240	-64.462	-3.077	-10.671	1.00	52.32	C
ATOM	11727	NH1	ARG B 240	-64.342	-4.141	-9.875	1.00	53.16	N
ATOM	11730	NH2	ARG B 240	-65.677	-2.577	-10.899	1.00	53.18	N
ATOM	11733	C	ARG B 240	-60.113	-2.914	-15.464	1.00	39.91	C
ATOM	11734	O	ARG B 240	-59.490	-2.000	-15.988	1.00	39.86	O
ATOM	11736	N	TRP B 241	-59.541	-4.002	-14.962	1.00	39.14	N
ATOM	11737	CA	TRP B 241	-58.105	-4.217	-15.004	1.00	38.37	C
ATOM	11739	CB	TRP B 241	-57.773	-5.635	-14.521	1.00	38.17	C
ATOM	11742	CG	TRP B 241	-56.373	-6.005	-14.773	1.00	36.58	C
ATOM	11743	CD1	TRP B 241	-55.290	-5.668	-14.022	1.00	36.32	C
ATOM	11745	NE1	TRP B 241	-54.152	-6.176	-14.587	1.00	35.63	N
ATOM	11747	CE2	TRP B 241	-54.494	-6.846	-15.731	1.00	34.31	C
ATOM	11748	CD2	TRP B 241	-55.883	-6.757	-15.874	1.00	33.99	C
ATOM	11749	CE3	TRP B 241	-56.484	-7.357	-16.977	1.00	33.28	C
ATOM	11751	CZ3	TRP B 241	-55.691	-8.034	-17.883	1.00	32.79	C
ATOM	11753	CH2	TRP B 241	-54.313	-8.101	-17.720	1.00	33.06	C
ATOM	11755	CZ2	TRP B 241	-53.697	-7.515	-16.648	1.00	34.16	C
ATOM	11757	C	TRP B 241	-57.561	-3.999	-16.409	1.00	38.25	C
ATOM	11758	O	TRP B 241	-56.607	-3.252	-16.603	1.00	38.13	O
ATOM	11760	N	TRP B 242	-58.186	-4.644	-17.386	1.00	38.21	N
ATOM	11761	CA	TRP B 242	-57.695	-4.628	-18.765	1.00	38.27	C
ATOM	11763	CB	TRP B 242	-58.479	-5.640	-19.609	1.00	38.31	C
ATOM	11766	CG	TRP B 242	-57.948	-5.872	-20.990	1.00	38.61	C
ATOM	11767	CD1	TRP B 242	-58.642	-5.754	-22.156	1.00	39.11	C
ATOM	11769	NE1	TRP B 242	-57.828	-6.046	-23.224	1.00	39.31	N
ATOM	11771	CE2	TRP B 242	-56.582	-6.363	-22.759	1.00	38.95	C

ATOM	11772	CD2 TRP B 242	-56.617	-6.263	-21.356	1.00	38.66	C
ATOM	11773	CE3 TRP B 242	-55.454	-6.528	-20.637	1.00	38.98	C
ATOM	11775	CZ3 TRP B 242	-54.314	-6.886	-21.326	1.00	39.20	C
ATOM	11777	CH2 TRP B 242	-54.313	-6.985	-22.716	1.00	39.12	C
ATOM	11779	CZ2 TRP B 242	-55.434	-6.724	-23.451	1.00	39.22	C
ATOM	11781	C TRP B 242	-57.739	-3.241	-19.400	1.00	38.27	C
ATOM	11782	O TRP B 242	-56.814	-2.873	-20.108	1.00	37.81	O
ATOM	11784	N ARG B 243	-58.809	-2.489	-19.140	1.00	38.79	N
ATOM	11785	CA ARG B 243	-58.936	-1.097	-19.607	1.00	39.37	C
ATOM	11787	CB ARG B 243	-60.351	-.552	-19.346	1.00	39.66	C
ATOM	11790	CG ARG B 243	-61.324	-.792	-20.520	1.00	41.96	C
ATOM	11793	CD ARG B 243	-62.687	-1.338	-20.085	1.00	44.82	C
ATOM	11796	NE ARG B 243	-63.461	-.380	-19.296	1.00	47.15	N
ATOM	11798	CZ ARG B 243	-64.389	-.704	-18.386	1.00	49.26	C
ATOM	11799	NH1 ARG B 243	-64.685	-1.976	-18.109	1.00	49.18	N
ATOM	11802	NH2 ARG B 243	-65.028	.262	-17.729	1.00	50.40	N
ATOM	11805	C ARG B 243	-57.890	-.183	-18.980	1.00	39.28	C
ATOM	11806	O ARG B 243	-57.246	.577	-19.681	1.00	39.40	O
ATOM	11808	N ARG B 244	-57.724	-.282	-17.665	1.00	39.52	N
ATOM	11809	CA ARG B 244	-56.734	.489	-16.888	1.00	39.61	C
ATOM	11811	CB ARG B 244	-56.774	.041	-15.409	1.00	40.25	C
ATOM	11814	CG ARG B 244	-56.115	.971	-14.383	1.00	42.01	C
ATOM	11817	CD ARG B 244	-57.029	2.154	-14.017	1.00	44.39	C
ATOM	11820	NE ARG B 244	-56.391	3.115	-13.103	1.00	46.36	N
ATOM	11822	CZ ARG B 244	-55.434	3.984	-13.441	1.00	48.06	C
ATOM	11823	NH1 ARG B 244	-54.951	4.043	-14.686	1.00	48.90	N
ATOM	11826	NH2 ARG B 244	-54.947	4.805	-12.519	1.00	48.82	N
ATOM	11829	C ARG B 244	-55.323	.322	-17.437	1.00	38.83	C
ATOM	11830	O ARG B 244	-54.594	1.300	-17.567	1.00	38.39	O
ATOM	11832	N VAL B 245	-54.953	-.922	-17.738	1.00	38.52	N
ATOM	11833	CA VAL B 245	-53.694	-1.241	-18.437	1.00	38.43	C
ATOM	11835	CB VAL B 245	-53.491	-2.760	-18.614	1.00	38.31	C
ATOM	11837	CG1 VAL B 245	-53.247	-3.426	-17.283	1.00	37.77	C
ATOM	11841	CG2 VAL B 245	-52.344	-3.036	-19.561	1.00	37.66	C
ATOM	11845	C VAL B 245	-53.650	-.638	-19.832	1.00	38.68	C
ATOM	11846	O VAL B 245	-52.615	-.153	-20.247	1.00	38.57	O
ATOM	11848	N GLY B 246	-54.765	-.724	-20.559	1.00	39.26	N
ATOM	11849	CA GLY B 246	-54.974	-.012	-21.834	1.00	39.62	C
ATOM	11852	C GLY B 246	-53.962	-.262	-22.942	1.00	40.08	C
ATOM	11853	O GLY B 246	-53.617	.650	-23.700	1.00	39.89	O
ATOM	11855	N LEU B 247	-53.502	-1.500	-23.062	1.00	40.75	N
ATOM	11856	CA LEU B 247	-52.347	-1.778	-23.905	1.00	41.34	C
ATOM	11858	CB LEU B 247	-51.655	-3.065	-23.459	1.00	41.15	C
ATOM	11861	CG LEU B 247	-50.132	-3.026	-23.298	1.00	40.49	C
ATOM	11863	CD1 LEU B 247	-49.622	-1.769	-22.605	1.00	39.50	C
ATOM	11867	CD2 LEU B 247	-49.695	-4.242	-22.526	1.00	39.98	C
ATOM	11871	C LEU B 247	-52.744	-1.836	-25.370	1.00	42.50	C
ATOM	11872	O LEU B 247	-52.005	-1.358	-26.225	1.00	42.33	O
ATOM	11874	N ALA B 248	-53.925	-2.390	-25.654	1.00	44.14	N
ATOM	11875	CA ALA B 248	-54.439	-2.494	-27.041	1.00	45.16	C
ATOM	11877	CB ALA B 248	-55.705	-3.347	-27.087	1.00	45.15	C
ATOM	11881	C ALA B 248	-54.692	-1.140	-27.739	1.00	45.98	C
ATOM	11882	O ALA B 248	-54.604	-1.053	-28.971	1.00	46.37	O
ATOM	11884	N THR B 249	-55.004	-.095	-26.975	1.00	46.64	N
ATOM	11885	CA THR B 249	-55.142	1.234	-27.570	1.00	47.25	C
ATOM	11887	CB THR B 249	-55.905	2.253	-26.655	1.00	47.44	C
ATOM	11889	OG1 THR B 249	-55.001	2.845	-25.706	1.00	47.52	O
ATOM	11891	CG2 THR B 249	-57.104	1.592	-25.928	1.00	47.44	C
ATOM	11895	C THR B 249	-53.756	1.793	-27.931	1.00	47.59	C
ATOM	11896	O THR B 249	-53.553	2.279	-29.049	1.00	48.07	O

ATOM	11898	N	LYS B 250	-52.808	1.710	-26.995	1.00	47.67	N
ATOM	11899	CA	LYS B 250	-51.469	2.299	-27.185	1.00	47.68	C
ATOM	11901	CB	LYS B 250	-50.793	2.553	-25.833	1.00	47.68	C
ATOM	11904	CG	LYS B 250	-51.428	3.673	-24.999	1.00	47.77	C
ATOM	11907	CD	LYS B 250	-51.142	5.086	-25.552	1.00	47.59	C
ATOM	11910	CE	LYS B 250	-49.676	5.496	-25.441	1.00	46.75	C
ATOM	11913	NZ	LYS B 250	-49.479	6.878	-25.948	1.00	46.49	N
ATOM	11917	C	LYS B 250	-50.541	1.465	-28.082	1.00	47.72	C
ATOM	11918	O	LYS B 250	-49.591	1.988	-28.663	1.00	47.54	O
ATOM	11920	N	LEU B 251	-50.804	.169	-28.184	1.00	47.82	N
ATOM	11921	CA	LEU B 251	-50.072	-.675	-29.118	1.00	47.96	C
ATOM	11923	CB	LEU B 251	-49.584	-1.970	-28.447	1.00	48.04	C
ATOM	11926	CG	LEU B 251	-48.109	-2.064	-28.033	1.00	47.39	C
ATOM	11928	CD1	LEU B 251	-47.659	-.881	-27.186	1.00	46.47	C
ATOM	11932	CD2	LEU B 251	-47.894	-3.376	-27.304	1.00	46.78	C
ATOM	11936	C	LEU B 251	-50.985	-.964	-30.306	1.00	48.19	C
ATOM	11937	O	LEU B 251	-51.975	-1.693	-30.194	1.00	48.38	O
ATOM	11939	N	HIS B 252	-50.627	-.396	-31.449	1.00	48.42	N
ATOM	11940	CA	HIS B 252	-51.521	-.329	-32.601	1.00	48.62	C
ATOM	11942	CB	HIS B 252	-51.110	.850	-33.492	1.00	48.99	C
ATOM	11945	CG	HIS B 252	-51.000	2.145	-32.742	1.00	50.56	C
ATOM	11946	ND1	HIS B 252	-52.095	2.935	-32.453	1.00	52.11	N
ATOM	11948	CE1	HIS B 252	-51.704	3.993	-31.766	1.00	52.68	C
ATOM	11950	NE2	HIS B 252	-50.397	3.912	-31.585	1.00	52.86	N
ATOM	11952	CD2	HIS B 252	-49.932	2.764	-32.181	1.00	51.68	C
ATOM	11954	C	HIS B 252	-51.582	-1.636	-33.383	1.00	48.14	C
ATOM	11955	O	HIS B 252	-52.615	-1.975	-33.935	1.00	48.02	O
ATOM	11957	N	PHE B 253	-50.481	-2.372	-33.405	1.00	48.10	N
ATOM	11958	CA	PHE B 253	-50.422	-3.702	-34.044	1.00	48.18	C
ATOM	11960	CB	PHE B 253	-48.965	-4.070	-34.375	1.00	48.16	C
ATOM	11963	CG	PHE B 253	-48.119	-4.311	-33.163	1.00	47.85	C
ATOM	11964	CD1	PHE B 253	-47.900	-5.596	-32.702	1.00	48.56	C
ATOM	11966	CE1	PHE B 253	-47.140	-5.812	-31.576	1.00	48.98	C
ATOM	11968	CZ	PHE B 253	-46.602	-4.729	-30.895	1.00	48.50	C
ATOM	11970	CE2	PHE B 253	-46.817	-3.454	-31.349	1.00	47.40	C
ATOM	11972	CD2	PHE B 253	-47.569	-3.249	-32.465	1.00	47.22	C
ATOM	11974	C	PHE B 253	-51.038	-4.819	-33.185	1.00	48.30	C
ATOM	11975	O	PHE B 253	-51.344	-5.903	-33.691	1.00	47.63	O
ATOM	11977	N	ALA B 254	-51.204	-4.541	-31.888	1.00	48.86	N
ATOM	11978	CA	ALA B 254	-51.653	-5.531	-30.903	1.00	49.13	C
ATOM	11980	CB	ALA B 254	-51.290	-5.071	-29.494	1.00	49.02	C
ATOM	11984	C	ALA B 254	-53.153	-5.806	-30.982	1.00	49.46	C
ATOM	11985	O	ALA B 254	-53.967	-4.870	-31.011	1.00	49.50	O
ATOM	11987	N	ARG B 255	-53.501	-7.097	-31.006	1.00	49.75	N
ATOM	11988	CA	ARG B 255	-54.896	-7.549	-30.918	1.00	49.89	C
ATOM	11990	CB	ARG B 255	-55.028	-9.051	-31.245	1.00	49.94	C
ATOM	11993	CG	ARG B 255	-54.839	-9.439	-32.724	1.00	49.93	C
ATOM	11996	CD	ARG B 255	-54.709	-10.967	-32.904	1.00	49.79	C
ATOM	11999	NE	ARG B 255	-53.527	-11.505	-32.219	1.00	50.11	N
ATOM	12001	CZ	ARG B 255	-53.276	-12.801	-32.010	1.00	50.54	C
ATOM	12002	NH1	ARG B 255	-54.117	-13.736	-32.439	1.00	51.36	N
ATOM	12005	NH2	ARG B 255	-52.175	-13.172	-31.360	1.00	50.21	N
ATOM	12008	C	ARG B 255	-55.425	-7.308	-29.511	1.00	49.83	C
ATOM	12009	O	ARG B 255	-54.666	-6.975	-28.599	1.00	49.95	O
ATOM	12011	N	ASP B 256	-56.736	-7.463	-29.350	1.00	49.72	N
ATOM	12012	CA	ASP B 256	-57.360	-7.483	-28.034	1.00	49.55	C
ATOM	12014	CB	ASP B 256	-58.183	-6.222	-27.812	1.00	49.52	C
ATOM	12017	CG	ASP B 256	-59.134	-6.362	-26.654	1.00	50.30	C
ATOM	12018	OD1	ASP B 256	-60.310	-5.944	-26.780	1.00	50.45	O
ATOM	12019	OD2	ASP B 256	-58.702	-6.926	-25.622	1.00	51.68	O

ATOM	12020	C	ASP B 256	-58.254	-8.721	-27.931	1.00	49.15	C
ATOM	12021	O	ASP B 256	-59.159	-8.901	-28.737	1.00	49.28	O
ATOM	12023	N	ARG B 257	-57.998	-9.565	-26.937	1.00	48.61	N
ATOM	12024	CA	ARG B 257	-58.743	-10.809	-26.765	1.00	48.20	C
ATOM	12026	CB	ARG B 257	-57.949	-11.984	-27.373	1.00	48.41	C
ATOM	12029	CG	ARG B 257	-57.505	-11.799	-28.821	1.00	49.37	C
ATOM	12032	CD	ARG B 257	-58.693	-11.636	-29.784	1.00	50.80	C
ATOM	12035	NE	ARG B 257	-58.908	-12.814	-30.627	1.00	52.11	N
ATOM	12037	CZ	ARG B 257	-58.678	-12.876	-31.944	1.00	52.99	C
ATOM	12038	NH1	ARG B 257	-58.214	-11.820	-32.623	1.00	52.37	N
ATOM	12041	NH2	ARG B 257	-58.918	-14.016	-32.596	1.00	52.97	N
ATOM	12044	C	ARG B 257	-59.032	-11.089	-25.280	1.00	47.32	C
ATOM	12045	O	ARG B 257	-58.579	-12.101	-24.735	1.00	47.46	O
ATOM	12047	N	LEU B 258	-59.770	-10.201	-24.616	1.00	45.94	N
ATOM	12048	CA	LEU B 258	-60.108	-10.430	-23.205	1.00	44.67	C
ATOM	12050	CB	LEU B 258	-60.332	-9.130	-22.438	1.00	44.61	C
ATOM	12053	CG	LEU B 258	-60.422	-9.323	-20.920	1.00	43.83	C
ATOM	12055	CD1	LEU B 258	-59.045	-9.552	-20.339	1.00	43.48	C
ATOM	12059	CD2	LEU B 258	-61.073	-8.139	-20.262	1.00	43.85	C
ATOM	12063	C	LEU B 258	-61.347	-11.286	-23.090	1.00	43.49	C
ATOM	12064	O	LEU B 258	-61.375	-12.220	-22.295	1.00	43.84	O
ATOM	12066	N	ILE B 259	-62.364	-10.968	-23.883	1.00	41.92	N
ATOM	12067	CA	ILE B 259	-63.617	-11.711	-23.841	1.00	40.80	C
ATOM	12069	CB	ILE B 259	-64.718	-11.132	-24.783	1.00	41.05	C
ATOM	12071	CG1	ILE B 259	-64.852	-9.598	-24.653	1.00	41.61	C
ATOM	12074	CD1	ILE B 259	-65.200	-8.874	-25.985	1.00	42.16	C
ATOM	12078	CG2	ILE B 259	-66.065	-11.816	-24.496	1.00	40.46	C
ATOM	12082	C	ILE B 259	-63.319	-13.151	-24.250	1.00	39.50	C
ATOM	12083	O	ILE B 259	-63.782	-14.082	-23.598	1.00	39.30	O
ATOM	12085	N	GLU B 260	-62.524	-13.318	-25.316	1.00	37.99	N
ATOM	12086	CA	GLU B 260	-62.120	-14.652	-25.812	1.00	36.60	C
ATOM	12088	CB	GLU B 260	-61.329	-14.574	-27.142	1.00	36.75	C
ATOM	12091	CG	GLU B 260	-62.167	-14.319	-28.425	1.00	37.79	C
ATOM	12094	CD	GLU B 260	-62.134	-12.851	-28.905	1.00	39.47	C
ATOM	12095	OE1	GLU B 260	-62.207	-11.923	-28.055	1.00	39.97	O
ATOM	12096	OE2	GLU B 260	-62.033	-12.628	-30.139	1.00	40.31	O
ATOM	12097	C	GLU B 260	-61.282	-15.381	-24.770	1.00	34.83	C
ATOM	12098	O	GLU B 260	-61.465	-16.568	-24.544	1.00	34.69	O
ATOM	12100	N	SER B 261	-60.365	-14.669	-24.131	1.00	32.89	N
ATOM	12101	CA	SER B 261	-59.508	-15.292	-23.139	1.00	31.48	C
ATOM	12103	CB	SER B 261	-58.341	-14.376	-22.768	1.00	31.55	C
ATOM	12106	OG	SER B 261	-57.298	-14.493	-23.728	1.00	31.75	O
ATOM	12108	C	SER B 261	-60.294	-15.716	-21.900	1.00	30.04	C
ATOM	12109	O	SER B 261	-59.921	-16.674	-21.208	1.00	29.70	O
ATOM	12111	N	PHE B 262	-61.389	-15.021	-21.618	1.00	28.25	N
ATOM	12112	CA	PHE B 262	-62.203	-15.404	-20.480	1.00	26.75	C
ATOM	12114	CB	PHE B 262	-63.111	-14.278	-20.009	1.00	26.36	C
ATOM	12117	CG	PHE B 262	-63.532	-14.453	-18.603	1.00	24.33	C
ATOM	12118	CD1	PHE B 262	-62.697	-14.092	-17.584	1.00	22.66	C
ATOM	12120	CE1	PHE B 262	-63.047	-14.283	-16.301	1.00	22.65	C
ATOM	12122	CZ	PHE B 262	-64.245	-14.873	-16.007	1.00	23.93	C
ATOM	12124	CE2	PHE B 262	-65.076	-15.267	-17.010	1.00	24.08	C
ATOM	12126	CD2	PHE B 262	-64.711	-15.067	-18.304	1.00	23.96	C
ATOM	12128	C	PHE B 262	-63.017	-16.660	-20.751	1.00	26.02	C
ATOM	12129	O	PHE B 262	-63.074	-17.546	-19.909	1.00	26.32	O
ATOM	12131	N	TYR B 263	-63.652	-16.732	-21.915	1.00	24.98	N
ATOM	12132	CA	TYR B 263	-64.299	-17.969	-22.396	1.00	24.14	C
ATOM	12134	CB	TYR B 263	-64.701	-17.756	-23.853	1.00	24.00	C
ATOM	12137	CG	TYR B 263	-65.080	-18.961	-24.667	1.00	24.32	C
ATOM	12138	CD1	TYR B 263	-66.302	-19.604	-24.489	1.00	25.03	C

ATOM	12140	CE1 TYR B 263	-66.666	-20.694	-25.285	1.00	25.81	C
ATOM	12142	CZ TYR B 263	-65.798	-21.132	-26.287	1.00	27.06	C
ATOM	12143	OH TYR B 263	-66.106	-22.191	-27.123	1.00	27.14	O
ATOM	12145	CE2 TYR B 263	-64.590	-20.489	-26.479	1.00	26.80	C
ATOM	12147	CD2 TYR B 263	-64.250	-19.404	-25.680	1.00	25.65	C
ATOM	12149	C TYR B 263	-63.351	-19.162	-22.256	1.00	23.35	C
ATOM	12150	O TYR B 263	-63.712	-20.219	-21.740	1.00	22.80	O
ATOM	12152	N TRP B 264	-62.119	-18.955	-22.702	1.00	22.64	N
ATOM	12153	CA TRP B 264	-61.047	-19.918	-22.523	1.00	21.95	C
ATOM	12155	CB TRP B 264	-59.737	-19.334	-23.061	1.00	21.84	C
ATOM	12158	CG TRP B 264	-58.603	-20.278	-22.964	1.00	22.90	C
ATOM	12159	CD1 TRP B 264	-57.700	-20.387	-21.943	1.00	23.98	C
ATOM	12161	NE1 TRP B 264	-56.800	-21.386	-22.214	1.00	23.95	N
ATOM	12163	CE2 TRP B 264	-57.124	-21.949	-23.417	1.00	23.51	C
ATOM	12164	CD2 TRP B 264	-58.252	-21.273	-23.915	1.00	23.40	C
ATOM	12165	CE3 TRP B 264	-58.776	-21.652	-25.147	1.00	23.68	C
ATOM	12167	CZ3 TRP B 264	-58.178	-22.669	-25.826	1.00	24.44	C
ATOM	12169	CH2 TRP B 264	-57.060	-23.324	-25.310	1.00	24.97	C
ATOM	12171	CZ2 TRP B 264	-56.521	-22.978	-24.103	1.00	24.44	C
ATOM	12173	C TRP B 264	-60.897	-20.326	-21.050	1.00	21.11	C
ATOM	12174	O TRP B 264	-60.768	-21.511	-20.748	1.00	20.98	O
ATOM	12176	N ALA B 265	-60.916	-19.344	-20.147	1.00	20.17	N
ATOM	12177	CA ALA B 265	-60.774	-19.602	-18.712	1.00	19.36	C
ATOM	12179	CB ALA B 265	-60.621	-18.305	-17.944	1.00	19.30	C
ATOM	12183	C ALA B 265	-61.947	-20.397	-18.171	1.00	18.63	C
ATOM	12184	O ALA B 265	-61.763	-21.349	-17.434	1.00	18.64	O
ATOM	12186	N VAL B 266	-63.159	-20.034	-18.557	1.00	18.11	N
ATOM	12187	CA VAL B 266	-64.335	-20.758	-18.075	1.00	17.60	C
ATOM	12189	CB VAL B 266	-65.652	-20.221	-18.675	1.00	17.38	C
ATOM	12191	CG1 VAL B 266	-65.902	-18.803	-18.199	1.00	17.09	C
ATOM	12195	CG2 VAL B 266	-66.817	-21.120	-18.312	1.00	15.40	C
ATOM	12199	C VAL B 266	-64.221	-22.246	-18.376	1.00	17.67	C
ATOM	12200	O VAL B 266	-64.766	-23.058	-17.647	1.00	18.27	O
ATOM	12202	N GLY B 267	-63.516	-22.607	-19.444	1.00	17.43	N
ATOM	12203	CA GLY B 267	-63.316	-24.013	-19.789	1.00	16.98	C
ATOM	12206	C GLY B 267	-62.253	-24.683	-18.947	1.00	16.55	C
ATOM	12207	O GLY B 267	-62.360	-25.859	-18.609	1.00	16.37	O
ATOM	12209	N VAL B 268	-61.215	-23.939	-18.608	1.00	16.30	N
ATOM	12210	CA VAL B 268	-60.160	-24.500	-17.794	1.00	16.23	C
ATOM	12212	CB VAL B 268	-58.893	-23.651	-17.840	1.00	15.68	C
ATOM	12214	CG1 VAL B 268	-57.857	-24.194	-16.894	1.00	15.01	C
ATOM	12218	CG2 VAL B 268	-58.361	-23.674	-19.235	1.00	15.01	C
ATOM	12222	C VAL B 268	-60.654	-24.718	-16.374	1.00	16.86	C
ATOM	12223	O VAL B 268	-60.421	-25.784	-15.813	1.00	16.93	O
ATOM	12225	N ALA B 269	-61.362	-23.733	-15.817	1.00	17.64	N
ATOM	12226	CA ALA B 269	-61.911	-23.838	-14.461	1.00	18.46	C
ATOM	12228	CB ALA B 269	-60.958	-23.215	-13.472	1.00	18.15	C
ATOM	12232	C ALA B 269	-63.312	-23.213	-14.340	1.00	19.34	C
ATOM	12233	O ALA B 269	-63.448	-22.030	-14.047	1.00	19.78	O
ATOM	12235	N PHE B 270	-64.347	-24.028	-14.529	1.00	20.34	N
ATOM	12236	CA PHE B 270	-65.724	-23.535	-14.631	1.00	21.21	C
ATOM	12238	CB PHE B 270	-66.564	-24.483	-15.482	1.00	21.34	C
ATOM	12241	CG PHE B 270	-67.083	-25.667	-14.718	1.00	22.46	C
ATOM	12242	CD1 PHE B 270	-68.294	-25.596	-14.039	1.00	23.19	C
ATOM	12244	CE1 PHE B 270	-68.756	-26.672	-13.315	1.00	23.27	C
ATOM	12246	CZ PHE B 270	-68.003	-27.836	-13.258	1.00	23.35	C
ATOM	12248	CE2 PHE B 270	-66.794	-27.918	-13.932	1.00	22.73	C
ATOM	12250	CD2 PHE B 270	-66.340	-26.843	-14.646	1.00	22.88	C
ATOM	12252	C PHE B 270	-66.464	-23.366	-13.316	1.00	21.98	C
ATOM	12253	O PHE B 270	-67.408	-22.583	-13.268	1.00	21.93	O

ATOM	12255	N	GLU B 271	-66.102	-24.126	-12.275	1.00	23.01	N
ATOM	12256	CA	GLU B 271	-66.929	-24.142	-11.037	1.00	23.99	C
ATOM	12258	CB	GLU B 271	-66.679	-25.353	-10.113	1.00	24.25	C
ATOM	12261	CG	GLU B 271	-65.286	-25.877	-10.092	1.00	26.22	C
ATOM	12264	CD	GLU B 271	-65.013	-26.914	-11.183	1.00	29.01	C
ATOM	12265	OE1	GLU B 271	-65.563	-28.032	-11.090	1.00	31.36	O
ATOM	12266	OE2	GLU B 271	-64.232	-26.622	-12.121	1.00	30.98	O
ATOM	12267	C	GLU B 271	-66.795	-22.831	-10.281	1.00	23.90	C
ATOM	12268	O	GLU B 271	-65.702	-22.296	-10.181	1.00	24.29	O
ATOM	12270	N	PRO B 272	-67.913	-22.316	-9.745	1.00	24.09	N
ATOM	12271	CA	PRO B 272	-68.005	-20.885	-9.428	1.00	24.11	C
ATOM	12273	CB	PRO B 272	-69.380	-20.751	-8.751	1.00	24.03	C
ATOM	12276	CG	PRO B 272	-70.083	-22.069	-8.997	1.00	24.12	C
ATOM	12279	CD	PRO B 272	-69.013	-23.083	-9.132	1.00	23.99	C
ATOM	12282	C	PRO B 272	-66.900	-20.386	-8.501	1.00	24.28	C
ATOM	12283	O	PRO B 272	-66.421	-19.261	-8.676	1.00	24.40	O
ATOM	12284	N	GLN B 273	-66.480	-21.222	-7.547	1.00	24.20	N
ATOM	12285	CA	GLN B 273	-65.500	-20.809	-6.551	1.00	24.20	C
ATOM	12287	CB	GLN B 273	-65.437	-21.840	-5.426	1.00	24.42	C
ATOM	12290	CG	GLN B 273	-64.750	-23.176	-5.775	1.00	25.17	C
ATOM	12293	CD	GLN B 273	-65.717	-24.334	-6.007	1.00	26.13	C
ATOM	12294	OE1	GLN B 273	-66.857	-24.141	-6.453	1.00	27.46	O
ATOM	12295	NE2	GLN B 273	-65.251	-25.554	-5.721	1.00	25.05	N
ATOM	12298	C	GLN B 273	-64.086	-20.546	-7.088	1.00	24.27	C
ATOM	12299	O	GLN B 273	-63.221	-20.151	-6.329	1.00	24.32	O
ATOM	12301	N	TYR B 274	-63.842	-20.759	-8.380	1.00	24.59	N
ATOM	12302	CA	TYR B 274	-62.509	-20.551	-8.962	1.00	24.85	C
ATOM	12304	CB	TYR B 274	-62.129	-21.724	-9.880	1.00	24.62	C
ATOM	12307	CG	TYR B 274	-62.009	-23.054	-9.184	1.00	24.35	C
ATOM	12308	CD1	TYR B 274	-61.395	-23.169	-7.946	1.00	24.62	C
ATOM	12310	CE1	TYR B 274	-61.280	-24.391	-7.314	1.00	25.10	C
ATOM	12312	CZ	TYR B 274	-61.771	-25.525	-7.920	1.00	24.92	C
ATOM	12313	OH	TYR B 274	-61.661	-26.748	-7.299	1.00	24.73	O
ATOM	12315	CE2	TYR B 274	-62.366	-25.437	-9.150	1.00	24.76	C
ATOM	12317	CD2	TYR B 274	-62.477	-24.203	-9.779	1.00	24.70	C
ATOM	12319	C	TYR B 274	-62.409	-19.235	-9.745	1.00	25.27	C
ATOM	12320	O	TYR B 274	-61.956	-19.210	-10.895	1.00	25.42	O
ATOM	12322	N	SER B 275	-62.808	-18.134	-9.121	1.00	25.40	N
ATOM	12323	CA	SER B 275	-62.755	-16.849	-9.802	1.00	25.43	C
ATOM	12325	CB	SER B 275	-63.615	-15.814	-9.071	1.00	25.53	C
ATOM	12328	OG	SER B 275	-65.006	-16.068	-9.291	1.00	25.52	O
ATOM	12330	C	SER B 275	-61.300	-16.403	-9.948	1.00	25.36	C
ATOM	12331	O	SER B 275	-60.885	-15.903	-10.996	1.00	25.17	O
ATOM	12333	N	ASP B 276	-60.517	-16.617	-8.901	1.00	25.41	N
ATOM	12334	CA	ASP B 276	-59.095	-16.353	-8.979	1.00	25.42	C
ATOM	12336	CB	ASP B 276	-58.390	-16.794	-7.712	1.00	25.44	C
ATOM	12339	CG	ASP B 276	-58.728	-15.925	-6.558	1.00	26.32	C
ATOM	12340	OD1	ASP B 276	-59.302	-14.857	-6.818	1.00	27.50	O
ATOM	12341	OD2	ASP B 276	-58.442	-16.299	-5.398	1.00	30.16	O
ATOM	12342	C	ASP B 276	-58.495	-17.069	-10.157	1.00	25.20	C
ATOM	12343	O	ASP B 276	-57.745	-16.461	-10.929	1.00	25.66	O
ATOM	12345	N	CYS B 277	-58.814	-18.353	-10.305	1.00	24.60	N
ATOM	12346	CA	CYS B 277	-58.228	-19.110	-11.391	1.00	24.30	C
ATOM	12348	CB	CYS B 277	-58.684	-20.561	-11.392	1.00	24.28	C
ATOM	12351	SG	CYS B 277	-57.737	-21.574	-12.569	1.00	23.82	S
ATOM	12353	C	CYS B 277	-58.593	-18.448	-12.698	1.00	24.09	C
ATOM	12354	O	CYS B 277	-57.727	-18.013	-13.450	1.00	23.87	O
ATOM	12356	N	ARG B 278	-59.893	-18.329	-12.933	1.00	24.22	N
ATOM	12357	CA	ARG B 278	-60.400	-17.746	-14.171	1.00	23.98	C
ATOM	12359	CB	ARG B 278	-61.917	-17.556	-14.123	1.00	23.93	C

ATOM	12362	CG	ARG B 278	-62.662	-18.871	-14.283	1.00	24.07	C
ATOM	12365	CD	ARG B 278	-64.132	-18.670	-14.607	1.00	24.48	C
ATOM	12368	NE	ARG B 278	-64.874	-18.141	-13.472	1.00	23.75	N
ATOM	12370	CZ	ARG B 278	-65.254	-18.856	-12.419	1.00	23.72	C
ATOM	12371	NH1	ARG B 278	-64.973	-20.149	-12.326	1.00	22.86	N
ATOM	12374	NH2	ARG B 278	-65.927	-18.262	-11.443	1.00	25.08	N
ATOM	12377	C	ARG B 278	-59.701	-16.453	-14.491	1.00	23.56	C
ATOM	12378	O	ARG B 278	-59.296	-16.265	-15.628	1.00	23.87	O
ATOM	12380	N	ASN B 279	-59.512	-15.593	-13.495	1.00	23.24	N
ATOM	12381	CA	ASN B 279	-58.854	-14.301	-13.734	1.00	23.48	C
ATOM	12383	CB	ASN B 279	-59.055	-13.340	-12.574	1.00	24.05	C
ATOM	12386	CG	ASN B 279	-60.510	-13.017	-12.346	1.00	26.58	C
ATOM	12387	OD1	ASN B 279	-61.378	-13.310	-13.189	1.00	28.96	O
ATOM	12388	ND2	ASN B 279	-60.799	-12.416	-11.194	1.00	30.30	N
ATOM	12391	C	ASN B 279	-57.376	-14.404	-14.027	1.00	22.58	C
ATOM	12392	O	ASN B 279	-56.903	-13.749	-14.941	1.00	22.73	O
ATOM	12394	N	SER B 280	-56.648	-15.205	-13.251	1.00	21.62	N
ATOM	12395	CA	SER B 280	-55.243	-15.483	-13.550	1.00	20.73	C
ATOM	12397	CB	SER B 280	-54.696	-16.597	-12.659	1.00	20.85	C
ATOM	12400	OG	SER B 280	-53.805	-16.099	-11.683	1.00	21.55	O
ATOM	12402	C	SER B 280	-55.091	-15.896	-15.004	1.00	20.04	C
ATOM	12403	O	SER B 280	-54.293	-15.312	-15.738	1.00	19.87	O
ATOM	12405	N	VAL B 281	-55.875	-16.893	-15.413	1.00	19.24	N
ATOM	12406	CA	VAL B 281	-55.720	-17.491	-16.720	1.00	18.77	C
ATOM	12408	CB	VAL B 281	-56.482	-18.834	-16.815	1.00	18.74	C
ATOM	12410	CG1	VAL B 281	-56.467	-19.404	-18.251	1.00	18.34	C
ATOM	12414	CG2	VAL B 281	-55.863	-19.844	-15.848	1.00	18.13	C
ATOM	12418	C	VAL B 281	-56.142	-16.486	-17.787	1.00	18.89	C
ATOM	12419	O	VAL B 281	-55.483	-16.364	-18.832	1.00	18.75	O
ATOM	12421	N	ALA B 282	-57.210	-15.739	-17.505	1.00	18.85	N
ATOM	12422	CA	ALA B 282	-57.709	-14.723	-18.439	1.00	18.80	C
ATOM	12424	CB	ALA B 282	-59.051	-14.191	-17.991	1.00	18.32	C
ATOM	12428	C	ALA B 282	-56.718	-13.576	-18.615	1.00	18.94	C
ATOM	12429	O	ALA B 282	-56.582	-13.021	-19.700	1.00	19.16	O
ATOM	12431	N	LYS B 283	-56.022	-13.211	-17.553	1.00	19.24	N
ATOM	12432	CA	LYS B 283	-55.037	-12.145	-17.658	1.00	19.62	C
ATOM	12434	CB	LYS B 283	-54.649	-11.615	-16.276	1.00	19.68	C
ATOM	12437	CG	LYS B 283	-55.779	-10.923	-15.518	1.00	19.95	C
ATOM	12440	CD	LYS B 283	-55.265	-10.398	-14.196	1.00	21.72	C
ATOM	12443	CE	LYS B 283	-56.364	-10.236	-13.164	1.00	23.83	C
ATOM	12446	NZ	LYS B 283	-56.129	-9.003	-12.351	1.00	24.72	N
ATOM	12450	C	LYS B 283	-53.802	-12.643	-18.396	1.00	19.77	C
ATOM	12451	O	LYS B 283	-53.237	-11.929	-19.208	1.00	20.05	O
ATOM	12453	N	MET B 284	-53.387	-13.871	-18.117	1.00	19.88	N
ATOM	12454	CA	MET B 284	-52.136	-14.367	-18.652	1.00	19.90	C
ATOM	12456	CB	MET B 284	-51.719	-15.676	-17.960	1.00	20.08	C
ATOM	12459	CG	MET B 284	-51.125	-15.505	-16.560	1.00	20.62	C
ATOM	12462	SD	MET B 284	-49.782	-14.295	-16.485	1.00	23.69	S
ATOM	12463	CE	MET B 284	-50.600	-12.872	-15.754	1.00	23.82	C
ATOM	12467	C	MET B 284	-52.257	-14.572	-20.145	1.00	19.82	C
ATOM	12468	O	MET B 284	-51.326	-14.269	-20.889	1.00	20.21	O
ATOM	12470	N	PHE B 285	-53.401	-15.088	-20.579	1.00	19.72	N
ATOM	12471	CA	PHE B 285	-53.620	-15.428	-21.986	1.00	19.66	C
ATOM	12473	CB	PHE B 285	-54.863	-16.324	-22.114	1.00	20.09	C
ATOM	12476	CG	PHE B 285	-54.986	-17.077	-23.427	1.00	21.18	C
ATOM	12477	CD1	PHE B 285	-53.877	-17.362	-24.232	1.00	22.86	C
ATOM	12479	CE1	PHE B 285	-54.013	-18.070	-25.424	1.00	22.58	C
ATOM	12481	CZ	PHE B 285	-55.254	-18.518	-25.806	1.00	23.49	C
ATOM	12483	CE2	PHE B 285	-56.365	-18.254	-25.002	1.00	23.36	C
ATOM	12485	CD2	PHE B 285	-56.222	-17.551	-23.823	1.00	21.76	C

ATOM	12487	C	PHE B 285	-53.781	-14.148	-22.781	1.00	19.18	C
ATOM	12488	O	PHE B 285	-53.333	-14.046	-23.914	1.00	18.98	O
ATOM	12490	N	SER B 286	-54.410	-13.161	-22.164	1.00	18.92	N
ATOM	12491	CA	SER B 286	-54.488	-11.843	-22.753	1.00	19.15	C
ATOM	12493	CB	SER B 286	-55.292	-10.894	-21.863	1.00	19.38	C
ATOM	12496	OG	SER B 286	-56.684	-11.180	-21.957	1.00	20.24	O
ATOM	12498	C	SER B 286	-53.106	-11.267	-23.032	1.00	18.91	C
ATOM	12499	O	SER B 286	-52.894	-10.668	-24.087	1.00	19.48	O
ATOM	12501	N	PHE B 287	-52.172	-11.446	-22.102	1.00	18.48	N
ATOM	12502	CA	PHE B 287	-50.768	-11.065	-22.339	1.00	18.30	C
ATOM	12504	CB	PHE B 287	-49.958	-11.035	-21.044	1.00	18.27	C
ATOM	12507	CG	PHE B 287	-50.040	-9.740	-20.324	1.00	18.59	C
ATOM	12508	CD1	PHE B 287	-49.375	-8.624	-20.813	1.00	20.02	C
ATOM	12510	CE1	PHE B 287	-49.456	-7.410	-20.152	1.00	20.57	C
ATOM	12512	CZ	PHE B 287	-50.208	-7.308	-18.985	1.00	20.02	C
ATOM	12514	CE2	PHE B 287	-50.869	-8.417	-18.503	1.00	19.78	C
ATOM	12516	CD2	PHE B 287	-50.781	-9.623	-19.164	1.00	19.05	C
ATOM	12518	C	PHE B 287	-50.067	-11.978	-23.339	1.00	18.01	C
ATOM	12519	O	PHE B 287	-49.374	-11.498	-24.229	1.00	18.00	O
ATOM	12521	N	VAL B 288	-50.236	-13.286	-23.194	1.00	17.69	N
ATOM	12522	CA	VAL B 288	-49.720	-14.206	-24.197	1.00	17.66	C
ATOM	12524	CB	VAL B 288	-50.168	-15.649	-23.930	1.00	17.37	C
ATOM	12526	CG1	VAL B 288	-49.937	-16.529	-25.134	1.00	16.24	C
ATOM	12530	CG2	VAL B 288	-49.414	-16.185	-22.758	1.00	17.21	C
ATOM	12534	C	VAL B 288	-50.081	-13.761	-25.633	1.00	18.08	C
ATOM	12535	O	VAL B 288	-49.207	-13.710	-26.492	1.00	18.08	O
ATOM	12537	N	THR B 289	-51.327	-13.396	-25.895	1.00	18.42	N
ATOM	12538	CA	THR B 289	-51.683	-13.001	-27.256	1.00	19.33	C
ATOM	12540	CB	THR B 289	-53.171	-12.736	-27.401	1.00	19.35	C
ATOM	12542	OG1	THR B 289	-53.573	-11.895	-26.321	1.00	21.59	O
ATOM	12544	CG2	THR B 289	-53.973	-14.046	-27.350	1.00	19.09	C
ATOM	12548	C	THR B 289	-50.901	-11.766	-27.723	1.00	19.71	C
ATOM	12549	O	THR B 289	-50.469	-11.707	-28.891	1.00	19.88	O
ATOM	12551	N	ILE B 290	-50.689	-10.796	-26.829	1.00	19.94	N
ATOM	12552	CA	ILE B 290	-49.943	-9.585	-27.225	1.00	20.39	C
ATOM	12554	CB	ILE B 290	-50.069	-8.410	-26.221	1.00	20.50	C
ATOM	12556	CG1	ILE B 290	-51.510	-8.137	-25.823	1.00	21.25	C
ATOM	12559	CD1	ILE B 290	-51.650	-6.854	-25.002	1.00	21.58	C
ATOM	12563	CG2	ILE B 290	-49.546	-7.134	-26.853	1.00	20.39	C
ATOM	12567	C	ILE B 290	-48.434	-9.840	-27.438	1.00	20.45	C
ATOM	12568	O	ILE B 290	-47.853	-9.317	-28.393	1.00	20.76	O
ATOM	12570	N	ILE B 291	-47.808	-10.616	-26.546	1.00	20.21	N
ATOM	12571	CA	ILE B 291	-46.375	-10.869	-26.618	1.00	19.87	C
ATOM	12573	CB	ILE B 291	-45.840	-11.603	-25.407	1.00	19.79	C
ATOM	12575	CG1	ILE B 291	-46.205	-10.879	-24.102	1.00	19.71	C
ATOM	12578	CD1	ILE B 291	-45.326	-9.734	-23.748	1.00	19.77	C
ATOM	12582	CG2	ILE B 291	-44.336	-11.747	-25.516	1.00	18.91	C
ATOM	12586	C	ILE B 291	-46.134	-11.735	-27.825	1.00	20.44	C
ATOM	12587	O	ILE B 291	-45.233	-11.481	-28.600	1.00	20.35	O
ATOM	12589	N	ASP B 292	-46.965	-12.751	-28.005	1.00	21.33	N
ATOM	12590	CA	ASP B 292	-46.915	-13.563	-29.228	1.00	22.13	C
ATOM	12592	CB	ASP B 292	-48.120	-14.508	-29.316	1.00	22.40	C
ATOM	12595	CG	ASP B 292	-47.847	-15.741	-30.160	1.00	23.47	C
ATOM	12596	OD1	ASP B 292	-46.852	-15.765	-30.923	1.00	24.83	O
ATOM	12597	OD2	ASP B 292	-48.637	-16.707	-30.038	1.00	25.22	O
ATOM	12598	C	ASP B 292	-46.856	-12.716	-30.507	1.00	22.28	C
ATOM	12599	O	ASP B 292	-46.088	-13.033	-31.399	1.00	22.73	O
ATOM	12601	N	ASP B 293	-47.656	-11.658	-30.611	1.00	22.22	N
ATOM	12602	CA	ASP B 293	-47.630	-10.828	-31.822	1.00	22.43	C
ATOM	12604	CB	ASP B 293	-48.840	-9.894	-31.885	1.00	23.03	C

ATOM	12607	CG	ASP B 293	-50.157	-10.644	-31.959	1.00	24.69	C
ATOM	12608	OD1	ASP B 293	-50.113	-11.900	-32.025	1.00	28.66	O
ATOM	12609	OD2	ASP B 293	-51.231	-9.986	-31.930	1.00	23.57	O
ATOM	12610	C	ASP B 293	-46.354	-9.997	-31.921	1.00	21.94	C
ATOM	12611	O	ASP B 293	-45.898	-9.674	-33.019	1.00	22.31	O
ATOM	12613	N	ILE B 294	-45.783	-9.643	-30.780	1.00	21.05	N
ATOM	12614	CA	ILE B 294	-44.513	-8.941	-30.770	1.00	20.53	C
ATOM	12616	CB	ILE B 294	-44.166	-8.448	-29.344	1.00	20.42	C
ATOM	12618	CG1	ILE B 294	-45.150	-7.349	-28.945	1.00	19.66	C
ATOM	12621	CD1	ILE B 294	-45.057	-6.956	-27.536	1.00	18.81	C
ATOM	12625	CG2	ILE B 294	-42.726	-7.924	-29.253	1.00	19.99	C
ATOM	12629	C	ILE B 294	-43.416	-9.845	-31.343	1.00	20.44	C
ATOM	12630	O	ILE B 294	-42.650	-9.431	-32.205	1.00	20.17	O
ATOM	12632	N	TYR B 295	-43.355	-11.088	-30.886	1.00	20.65	N
ATOM	12633	CA	TYR B 295	-42.313	-12.008	-31.353	1.00	20.76	C
ATOM	12635	CB	TYR B 295	-42.145	-13.198	-30.412	1.00	20.29	C
ATOM	12638	CG	TYR B 295	-41.393	-12.926	-29.129	1.00	17.92	C
ATOM	12639	CD1	TYR B 295	-40.117	-13.409	-28.936	1.00	17.62	C
ATOM	12641	CE1	TYR B 295	-39.427	-13.195	-27.747	1.00	17.12	C
ATOM	12643	CZ	TYR B 295	-40.022	-12.493	-26.730	1.00	16.61	C
ATOM	12644	OH	TYR B 295	-39.357	-12.281	-25.528	1.00	14.46	O
ATOM	12646	CE2	TYR B 295	-41.296	-12.008	-26.911	1.00	16.82	C
ATOM	12648	CD2	TYR B 295	-41.972	-12.235	-28.102	1.00	16.39	C
ATOM	12650	C	TYR B 295	-42.619	-12.519	-32.749	1.00	21.65	C
ATOM	12651	O	TYR B 295	-41.698	-12.751	-33.540	1.00	22.13	O
ATOM	12653	N	ASP B 296	-43.905	-12.677	-33.049	1.00	22.47	N
ATOM	12654	CA	ASP B 296	-44.327	-13.316	-34.285	1.00	23.52	C
ATOM	12656	CB	ASP B 296	-45.761	-13.851	-34.185	1.00	23.99	C
ATOM	12659	CG	ASP B 296	-46.134	-14.744	-35.369	1.00	26.22	C
ATOM	12660	OD1	ASP B 296	-45.508	-15.832	-35.538	1.00	27.70	O
ATOM	12661	OD2	ASP B 296	-47.052	-14.348	-36.128	1.00	28.78	O
ATOM	12662	C	ASP B 296	-44.228	-12.392	-35.478	1.00	23.68	C
ATOM	12663	O	ASP B 296	-43.673	-12.784	-36.502	1.00	24.01	O
ATOM	12665	N	VAL B 297	-44.765	-11.178	-35.364	1.00	23.90	N
ATOM	12666	CA	VAL B 297	-44.801	-10.273	-36.524	1.00	23.99	C
ATOM	12668	CB	VAL B 297	-46.251	-10.003	-36.973	1.00	23.87	C
ATOM	12670	CG1	VAL B 297	-46.949	-11.310	-37.197	1.00	23.75	C
ATOM	12674	CG2	VAL B 297	-47.002	-9.136	-35.966	1.00	23.05	C
ATOM	12678	C	VAL B 297	-44.060	-8.940	-36.383	1.00	24.27	C
ATOM	12679	O	VAL B 297	-43.485	-8.465	-37.342	1.00	23.84	O
ATOM	12681	N	TYR B 298	-44.058	-8.347	-35.199	1.00	24.87	N
ATOM	12682	CA	TYR B 298	-43.760	-6.926	-35.078	1.00	25.67	C
ATOM	12684	CB	TYR B 298	-44.656	-6.294	-34.010	1.00	26.06	C
ATOM	12687	CG	TYR B 298	-44.520	-4.789	-33.932	1.00	27.95	C
ATOM	12688	CD1	TYR B 298	-45.097	-3.968	-34.901	1.00	29.74	C
ATOM	12690	CE1	TYR B 298	-44.982	-2.574	-34.843	1.00	29.93	C
ATOM	12692	CZ	TYR B 298	-44.282	-1.991	-33.806	1.00	29.96	C
ATOM	12693	OH	TYR B 298	-44.165	-.631	-33.746	1.00	30.01	O
ATOM	12695	CE2	TYR B 298	-43.697	-2.777	-32.829	1.00	30.39	C
ATOM	12697	CD2	TYR B 298	-43.811	-4.181	-32.899	1.00	29.96	C
ATOM	12699	C	TYR B 298	-42.318	-6.592	-34.755	1.00	25.76	C
ATOM	12700	O	TYR B 298	-41.696	-5.830	-35.470	1.00	26.54	O
ATOM	12702	N	GLY B 299	-41.798	-7.111	-33.653	1.00	25.89	N
ATOM	12703	CA	GLY B 299	-40.473	-6.706	-33.167	1.00	25.78	C
ATOM	12706	C	GLY B 299	-39.346	-7.341	-33.951	1.00	25.52	C
ATOM	12707	O	GLY B 299	-39.513	-8.410	-34.525	1.00	26.04	O
ATOM	12709	N	THR B 300	-38.196	-6.685	-33.990	1.00	25.32	N
ATOM	12710	CA	THR B 300	-37.063	-7.237	-34.716	1.00	25.16	C
ATOM	12712	CB	THR B 300	-36.088	-6.164	-35.295	1.00	25.08	C
ATOM	12714	OG1	THR B 300	-35.191	-5.710	-34.281	1.00	25.07	O

ATOM	12716	CG2 THR B 300	-36.833	-4.978	-35.875	1.00	25.08	C
ATOM	12720	C THR B 300	-36.324	-8.139	-33.771	1.00	25.13	C
ATOM	12721	O THR B 300	-36.335	-7.950	-32.570	1.00	25.28	O
ATOM	12723	N LEU B 301	-35.662	-9.129	-34.332	1.00	25.50	N
ATOM	12724	CA LEU B 301	-34.918	-10.093	-33.538	1.00	25.42	C
ATOM	12726	CB LEU B 301	-34.148	-11.020	-34.479	1.00	25.17	C
ATOM	12729	CG LEU B 301	-33.624	-12.338	-33.941	1.00	24.95	C
ATOM	12731	CD1 LEU B 301	-34.715	-13.233	-33.437	1.00	24.42	C
ATOM	12735	CD2 LEU B 301	-32.900	-12.998	-35.080	1.00	26.50	C
ATOM	12739	C LEU B 301	-33.995	-9.425	-32.485	1.00	25.54	C
ATOM	12740	O LEU B 301	-33.843	-9.959	-31.394	1.00	25.35	O
ATOM	12742	N ASP B 302	-33.417	-8.257	-32.782	1.00	25.68	N
ATOM	12743	CA ASP B 302	-32.598	-7.550	-31.774	1.00	26.02	C
ATOM	12745	CB ASP B 302	-31.800	-6.377	-32.378	1.00	26.37	C
ATOM	12748	CG ASP B 302	-30.563	-6.833	-33.167	1.00	28.27	C
ATOM	12749	OD1 ASP B 302	-29.508	-6.156	-33.048	1.00	29.05	O
ATOM	12750	OD2 ASP B 302	-30.652	-7.852	-33.910	1.00	31.04	O
ATOM	12751	C ASP B 302	-33.472	-7.037	-30.630	1.00	25.69	C
ATOM	12752	O ASP B 302	-33.110	-7.142	-29.447	1.00	25.90	O
ATOM	12754	N GLU B 303	-34.622	-6.473	-30.993	1.00	25.04	N
ATOM	12755	CA GLU B 303	-35.600	-6.029	-30.011	1.00	24.33	C
ATOM	12757	CB GLU B 303	-36.783	-5.346	-30.697	1.00	24.20	C
ATOM	12760	CG GLU B 303	-36.399	-4.109	-31.466	1.00	24.23	C
ATOM	12763	CD GLU B 303	-37.589	-3.382	-32.063	1.00	24.23	C
ATOM	12764	OE1 GLU B 303	-38.496	-4.065	-32.604	1.00	22.66	O
ATOM	12765	OE2 GLU B 303	-37.592	-2.121	-31.995	1.00	23.93	O
ATOM	12766	C GLU B 303	-36.103	-7.212	-29.193	1.00	23.87	C
ATOM	12767	O GLU B 303	-36.396	-7.058	-28.011	1.00	23.59	O
ATOM	12769	N LEU B 304	-36.204	-8.383	-29.825	1.00	23.39	N
ATOM	12770	CA LEU B 304	-36.766	-9.570	-29.166	1.00	23.23	C
ATOM	12772	CB LEU B 304	-37.180	-10.624	-30.208	1.00	22.83	C
ATOM	12775	CG LEU B 304	-38.381	-10.217	-31.074	1.00	22.06	C
ATOM	12777	CD1 LEU B 304	-38.822	-11.321	-32.037	1.00	20.11	C
ATOM	12781	CD2 LEU B 304	-39.551	-9.800	-30.163	1.00	21.94	C
ATOM	12785	C LEU B 304	-35.810	-10.147	-28.106	1.00	23.37	C
ATOM	12786	O LEU B 304	-36.248	-10.785	-27.141	1.00	23.22	O
ATOM	12788	N GLU B 305	-34.514	-9.892	-28.283	1.00	23.30	N
ATOM	12789	CA GLU B 305	-33.512	-10.308	-27.324	1.00	23.18	C
ATOM	12791	CB GLU B 305	-32.114	-10.176	-27.899	1.00	23.58	C
ATOM	12794	CG GLU B 305	-31.832	-11.142	-29.047	1.00	25.27	C
ATOM	12797	CD GLU B 305	-31.306	-12.488	-28.585	1.00	27.91	C
ATOM	12798	OE1 GLU B 305	-30.946	-13.303	-29.470	1.00	29.87	O
ATOM	12799	OE2 GLU B 305	-31.244	-12.728	-27.351	1.00	29.28	O
ATOM	12800	C GLU B 305	-33.633	-9.440	-26.104	1.00	22.52	C
ATOM	12801	O GLU B 305	-33.793	-9.956	-25.005	1.00	23.39	O
ATOM	12803	N LEU B 306	-33.576	-8.124	-26.274	1.00	21.51	N
ATOM	12804	CA LEU B 306	-33.715	-7.224	-25.114	1.00	20.70	C
ATOM	12806	CB LEU B 306	-33.899	-5.764	-25.544	1.00	20.59	C
ATOM	12809	CG LEU B 306	-32.769	-5.065	-26.293	1.00	20.12	C
ATOM	12811	CD1 LEU B 306	-32.992	-3.576	-26.203	1.00	19.89	C
ATOM	12815	CD2 LEU B 306	-31.402	-5.437	-25.739	1.00	19.81	C
ATOM	12819	C LEU B 306	-34.896	-7.645	-24.237	1.00	19.91	C
ATOM	12820	O LEU B 306	-34.771	-7.778	-23.031	1.00	19.53	O
ATOM	12822	N PHE B 307	-36.030	-7.883	-24.872	1.00	19.48	N
ATOM	12823	CA PHE B 307	-37.231	-8.267	-24.170	1.00	19.36	C
ATOM	12825	CB PHE B 307	-38.400	-8.406	-25.138	1.00	19.54	C
ATOM	12828	CG PHE B 307	-39.729	-8.296	-24.482	1.00	19.66	C
ATOM	12829	CD1 PHE B 307	-40.365	-7.079	-24.400	1.00	21.03	C
ATOM	12831	CE1 PHE B 307	-41.590	-6.968	-23.787	1.00	21.62	C
ATOM	12833	CZ PHE B 307	-42.182	-8.072	-23.248	1.00	20.25	C

ATOM	12835	CE2 PHE B 307	-41.553	-9.292	-23.329	1.00	20.07	C
ATOM	12837	CD2 PHE B 307	-40.338	-9.400	-23.939	1.00	19.62	C
ATOM	12839	C PHE B 307	-37.039	-9.577	-23.443	1.00	19.22	C
ATOM	12840	O PHE B 307	-37.354	-9.684	-22.260	1.00	19.50	O
ATOM	12842	N THR B 308	-36.532	-10.583	-24.143	1.00	18.96	N
ATOM	12843	CA THR B 308	-36.321	-11.874	-23.511	1.00	18.84	C
ATOM	12845	CB THR B 308	-35.774	-12.879	-24.490	1.00	18.62	C
ATOM	12847	OG1 THR B 308	-36.687	-12.988	-25.579	1.00	18.88	O
ATOM	12849	CG2 THR B 308	-35.626	-14.228	-23.833	1.00	18.74	C
ATOM	12853	C THR B 308	-35.384	-11.737	-22.317	1.00	19.04	C
ATOM	12854	O THR B 308	-35.686	-12.244	-21.233	1.00	18.69	O
ATOM	12856	N ASP B 309	-34.264	-11.031	-22.511	1.00	19.36	N
ATOM	12857	CA ASP B 309	-33.335	-10.751	-21.421	1.00	19.65	C
ATOM	12859	CB ASP B 309	-32.121	-9.983	-21.942	1.00	19.92	C
ATOM	12862	CG ASP B 309	-31.171	-9.536	-20.814	1.00	23.51	C
ATOM	12863	OD1 ASP B 309	-30.520	-10.412	-20.175	1.00	27.27	O
ATOM	12864	OD2 ASP B 309	-31.076	-8.299	-20.561	1.00	27.40	O
ATOM	12865	C ASP B 309	-34.048	-9.977	-20.295	1.00	19.08	C
ATOM	12866	O ASP B 309	-33.912	-10.308	-19.124	1.00	18.55	O
ATOM	12868	N ALA B 310	-34.836	-8.975	-20.665	1.00	18.93	N
ATOM	12869	CA ALA B 310	-35.546	-8.151	-19.688	1.00	19.12	C
ATOM	12871	CB ALA B 310	-36.346	-7.062	-20.380	1.00	18.79	C
ATOM	12875	C ALA B 310	-36.462	-8.965	-18.784	1.00	19.36	C
ATOM	12876	O ALA B 310	-36.560	-8.696	-17.577	1.00	19.58	O
ATOM	12878	N VAL B 311	-37.143	-9.942	-19.370	1.00	19.61	N
ATOM	12879	CA VAL B 311	-37.996	-10.845	-18.611	1.00	19.68	C
ATOM	12881	CB VAL B 311	-38.867	-11.682	-19.555	1.00	19.46	C
ATOM	12883	CG1 VAL B 311	-39.948	-10.825	-20.123	1.00	18.49	C
ATOM	12887	CG2 VAL B 311	-39.469	-12.863	-18.843	1.00	19.31	C
ATOM	12891	C VAL B 311	-37.161	-11.726	-17.673	1.00	20.46	C
ATOM	12892	O VAL B 311	-37.479	-11.843	-16.494	1.00	20.25	O
ATOM	12894	N GLU B 312	-36.085	-12.317	-18.189	1.00	21.67	N
ATOM	12895	CA GLU B 312	-35.198	-13.168	-17.384	1.00	22.65	C
ATOM	12897	CB GLU B 312	-34.056	-13.708	-18.239	1.00	23.03	C
ATOM	12900	CG GLU B 312	-34.454	-14.764	-19.261	1.00	24.88	C
ATOM	12903	CD GLU B 312	-33.335	-15.077	-20.267	1.00	27.34	C
ATOM	12904	OE1 GLU B 312	-33.485	-16.057	-21.028	1.00	29.07	O
ATOM	12905	OE2 GLU B 312	-32.316	-14.342	-20.313	1.00	28.63	O
ATOM	12906	C GLU B 312	-34.598	-12.443	-16.167	1.00	23.09	C
ATOM	12907	O GLU B 312	-34.506	-13.007	-15.082	1.00	23.04	O
ATOM	12909	N ARG B 313	-34.186	-11.199	-16.335	1.00	23.67	N
ATOM	12910	CA ARG B 313	-33.590	-10.486	-15.222	1.00	24.60	C
ATOM	12912	CB ARG B 313	-32.609	-9.436	-15.735	1.00	25.26	C
ATOM	12915	CG ARG B 313	-31.333	-10.030	-16.398	1.00	27.83	C
ATOM	12918	CD ARG B 313	-30.289	-8.943	-16.754	1.00	31.87	C
ATOM	12921	NE ARG B 313	-30.922	-7.793	-17.430	1.00	35.71	N
ATOM	12923	CZ ARG B 313	-31.341	-6.661	-16.837	1.00	38.42	C
ATOM	12924	NH1 ARG B 313	-31.189	-6.452	-15.522	1.00	39.15	N
ATOM	12927	NH2 ARG B 313	-31.919	-5.711	-17.576	1.00	39.25	N
ATOM	12930	C ARG B 313	-34.638	-9.879	-14.280	1.00	24.67	C
ATOM	12931	O ARG B 313	-34.359	-9.645	-13.117	1.00	24.09	O
ATOM	12933	N TRP B 314	-35.843	-9.635	-14.781	1.00	25.60	N
ATOM	12934	CA TRP B 314	-36.971	-9.180	-13.951	1.00	26.20	C
ATOM	12936	CB TRP B 314	-37.488	-10.339	-13.097	1.00	25.90	C
ATOM	12939	CG TRP B 314	-38.912	-10.184	-12.662	1.00	23.93	C
ATOM	12940	CD1 TRP B 314	-39.355	-9.854	-11.424	1.00	22.76	C
ATOM	12942	NE1 TRP B 314	-40.727	-9.804	-11.412	1.00	21.93	N
ATOM	12944	CE2 TRP B 314	-41.189	-10.107	-12.659	1.00	20.38	C
ATOM	12945	CD2 TRP B 314	-40.074	-10.352	-13.473	1.00	21.17	C
ATOM	12946	CE3 TRP B 314	-40.277	-10.684	-14.811	1.00	20.19	C

ATOM	12948	CZ3 TRP B 314	-41.549	-10.766	-15.280	1.00	19.79	C
ATOM	12950	CH2 TRP B 314	-42.640	-10.512	-14.445	1.00	21.00	C
ATOM	12952	CZ2 TRP B 314	-42.474	-10.179	-13.130	1.00	20.60	C
ATOM	12954	C TRP B 314	-36.620	-7.966	-13.076	1.00	27.57	C
ATOM	12955	O TRP B 314	-36.890	-7.935	-11.870	1.00	27.56	O
ATOM	12957	N ASP B 315	-36.039	-6.959	-13.718	1.00	29.12	N
ATOM	12958	CA ASP B 315	-35.452	-5.822	-13.027	1.00	30.46	C
ATOM	12960	CB ASP B 315	-33.935	-5.796	-13.315	1.00	30.84	C
ATOM	12963	CG ASP B 315	-33.227	-4.551	-12.772	1.00	32.75	C
ATOM	12964	OD1 ASP B 315	-33.747	-3.904	-11.830	1.00	35.67	O
ATOM	12965	OD2 ASP B 315	-32.128	-4.222	-13.295	1.00	34.91	O
ATOM	12966	C ASP B 315	-36.157	-4.592	-13.557	1.00	31.10	C
ATOM	12967	O ASP B 315	-35.888	-4.157	-14.674	1.00	31.35	O
ATOM	12969	N VAL B 316	-37.088	-4.044	-12.781	1.00	32.06	N
ATOM	12970	CA VAL B 316	-37.864	-2.896	-13.262	1.00	32.65	C
ATOM	12972	CB VAL B 316	-39.016	-2.510	-12.334	1.00	32.45	C
ATOM	12974	CG1 VAL B 316	-39.911	-1.499	-13.012	1.00	32.15	C
ATOM	12978	CG2 VAL B 316	-38.490	-1.960	-11.040	1.00	32.91	C
ATOM	12982	C VAL B 316	-36.964	-1.692	-13.442	1.00	33.43	C
ATOM	12983	O VAL B 316	-37.228	-.854	-14.294	1.00	33.74	O
ATOM	12985	N ASN B 317	-35.881	-1.628	-12.663	1.00	34.29	N
ATOM	12986	CA ASN B 317	-34.937	-.510	-12.730	1.00	34.63	C
ATOM	12988	CB ASN B 317	-33.995	-.526	-11.505	1.00	34.50	C
ATOM	12991	CG ASN B 317	-34.692	-.084	-10.201	1.00	34.08	C
ATOM	12992	OD1 ASN B 317	-35.198	1.042	-10.104	1.00	33.72	O
ATOM	12993	ND2 ASN B 317	-34.690	-.961	-9.193	1.00	31.22	N
ATOM	12996	C ASN B 317	-34.124	-.481	-14.029	1.00	35.21	C
ATOM	12997	O ASN B 317	-33.196	.296	-14.127	1.00	35.52	O
ATOM	12999	N ALA B 318	-34.470	-1.318	-15.013	1.00	35.89	N
ATOM	13000	CA ALA B 318	-33.757	-1.382	-16.299	1.00	36.48	C
ATOM	13002	CB ALA B 318	-32.727	-2.508	-16.253	1.00	36.49	C
ATOM	13006	C ALA B 318	-34.695	-1.525	-17.534	1.00	37.02	C
ATOM	13007	O ALA B 318	-34.297	-1.984	-18.616	1.00	36.52	O
ATOM	13009	N ILE B 319	-35.952	-1.137	-17.339	1.00	37.76	N
ATOM	13010	CA ILE B 319	-36.848	-.733	-18.424	1.00	38.08	C
ATOM	13012	CB ILE B 319	-37.926	.247	-17.893	1.00	38.11	C
ATOM	13014	CG1 ILE B 319	-39.110	-.495	-17.293	1.00	38.08	C
ATOM	13017	CD1 ILE B 319	-40.068	.451	-16.596	1.00	38.27	C
ATOM	13021	CG2 ILE B 319	-38.408	1.189	-18.996	1.00	37.84	C
ATOM	13025	C ILE B 319	-36.134	.053	-19.516	1.00	38.29	C
ATOM	13026	O ILE B 319	-36.192	-.324	-20.677	1.00	38.68	O
ATOM	13028	N ASN B 320	-35.460	1.141	-19.125	1.00	38.36	N
ATOM	13029	CA ASN B 320	-34.977	2.168	-20.074	1.00	38.17	C
ATOM	13031	CB ASN B 320	-34.278	3.331	-19.331	1.00	38.16	C
ATOM	13034	CG ASN B 320	-35.258	4.207	-18.514	1.00	37.87	C
ATOM	13035	OD1 ASN B 320	-36.366	4.542	-18.955	1.00	36.50	O
ATOM	13036	ND2 ASN B 320	-34.827	4.589	-17.322	1.00	38.12	N
ATOM	13039	C ASN B 320	-34.082	1.643	-21.206	1.00	37.78	C
ATOM	13040	O ASN B 320	-33.856	2.332	-22.181	1.00	37.76	O
ATOM	13042	N ASP B 321	-33.600	.416	-21.083	1.00	37.46	N
ATOM	13043	CA ASP B 321	-32.852	-.224	-22.162	1.00	37.25	C
ATOM	13045	CB ASP B 321	-32.091	-1.467	-21.629	1.00	37.83	C
ATOM	13048	CG ASP B 321	-31.351	-1.207	-20.285	1.00	39.33	C
ATOM	13049	OD1 ASP B 321	-31.247	-.021	-19.871	1.00	40.72	O
ATOM	13050	OD2 ASP B 321	-30.887	-2.199	-19.648	1.00	39.62	O
ATOM	13051	C ASP B 321	-33.764	-.636	-23.339	1.00	35.99	C
ATOM	13052	O ASP B 321	-33.273	-.904	-24.436	1.00	35.87	O
ATOM	13054	N LEU B 322	-35.076	-.674	-23.111	1.00	34.66	N
ATOM	13055	CA LEU B 322	-36.024	-1.263	-24.061	1.00	33.77	C
ATOM	13057	CB LEU B 322	-37.180	-1.961	-23.327	1.00	33.65	C

ATOM	13060	CG	LEU	B	322	-36.994	-3.260	-22.544	1.00	32.41	C
ATOM	13062	CD1	LEU	B	322	-38.223	-3.482	-21.693	1.00	31.54	C
ATOM	13066	CD2	LEU	B	322	-36.769	-4.432	-23.460	1.00	30.73	C
ATOM	13070	C	LEU	B	322	-36.670	-.243	-24.982	1.00	33.32	C
ATOM	13071	O	LEU	B	322	-36.910	.893	-24.570	1.00	33.37	O
ATOM	13073	N	PRO	B	323	-37.005	-.667	-26.218	1.00	32.72	N
ATOM	13074	CA	PRO	B	323	-37.856	.068	-27.126	1.00	32.42	C
ATOM	13076	CB	PRO	B	323	-38.239	-.985	-28.158	1.00	32.28	C
ATOM	13079	CG	PRO	B	323	-37.098	-1.849	-28.231	1.00	32.46	C
ATOM	13082	CD	PRO	B	323	-36.498	-1.891	-26.857	1.00	32.79	C
ATOM	13085	C	PRO	B	323	-39.115	.576	-26.453	1.00	32.50	C
ATOM	13086	O	PRO	B	323	-39.655	-.083	-25.550	1.00	32.46	O
ATOM	13087	N	ASP	B	324	-39.590	1.725	-26.923	1.00	32.44	N
ATOM	13088	CA	ASP	B	324	-40.729	2.387	-26.318	1.00	32.44	C
ATOM	13090	CB	ASP	B	324	-41.059	3.678	-27.071	1.00	32.91	C
ATOM	13093	CG	ASP	B	324	-40.241	4.873	-26.576	1.00	34.17	C
ATOM	13094	OD1	ASP	B	324	-39.507	4.728	-25.562	1.00	35.83	O
ATOM	13095	OD2	ASP	B	324	-40.351	5.960	-27.196	1.00	35.42	O
ATOM	13096	C	ASP	B	324	-41.967	1.504	-26.186	1.00	31.78	C
ATOM	13097	O	ASP	B	324	-42.546	1.440	-25.102	1.00	31.77	O
ATOM	13099	N	TYR	B	325	-42.368	.821	-27.258	1.00	30.99	N
ATOM	13100	CA	TYR	B	325	-43.519	-.086	-27.167	1.00	30.49	C
ATOM	13102	CB	TYR	B	325	-43.941	-.640	-28.536	1.00	30.28	C
ATOM	13105	CG	TYR	B	325	-43.027	-1.688	-29.134	1.00	30.25	C
ATOM	13106	CD1	TYR	B	325	-41.917	-1.332	-29.893	1.00	29.81	C
ATOM	13108	CE1	TYR	B	325	-41.089	-2.300	-30.449	1.00	29.74	C
ATOM	13110	CZ	TYR	B	325	-41.371	-3.643	-30.258	1.00	30.22	C
ATOM	13111	OH	TYR	B	325	-40.565	-4.625	-30.805	1.00	31.28	O
ATOM	13113	CE2	TYR	B	325	-42.466	-4.015	-29.519	1.00	30.23	C
ATOM	13115	CD2	TYR	B	325	-43.288	-3.038	-28.963	1.00	30.49	C
ATOM	13117	C	TYR	B	325	-43.284	-1.221	-26.168	1.00	30.14	C
ATOM	13118	O	TYR	B	325	-44.219	-1.638	-25.483	1.00	30.19	O
ATOM	13120	N	MET	B	326	-42.048	-1.700	-26.060	1.00	29.74	N
ATOM	13121	CA	MET	B	326	-41.755	-2.820	-25.155	1.00	29.82	C
ATOM	13123	CB	MET	B	326	-40.461	-3.541	-25.549	1.00	29.69	C
ATOM	13126	CG	MET	B	326	-40.650	-4.510	-26.686	1.00	28.83	C
ATOM	13129	SD	MET	B	326	-39.123	-5.272	-27.204	1.00	26.99	S
ATOM	13130	CE	MET	B	326	-39.761	-6.499	-28.328	1.00	29.13	C
ATOM	13134	C	MET	B	326	-41.696	-2.402	-23.683	1.00	29.92	C
ATOM	13135	O	MET	B	326	-42.032	-3.204	-22.801	1.00	29.75	O
ATOM	13137	N	LYS	B	327	-41.251	-1.166	-23.429	1.00	29.68	N
ATOM	13138	CA	LYS	B	327	-41.318	-.578	-22.088	1.00	29.46	C
ATOM	13140	CB	LYS	B	327	-41.016	.921	-22.114	1.00	29.93	C
ATOM	13143	CG	LYS	B	327	-39.572	1.305	-21.951	1.00	31.63	C
ATOM	13146	CD	LYS	B	327	-39.410	2.803	-22.249	1.00	34.93	C
ATOM	13149	CE	LYS	B	327	-37.986	3.291	-22.021	1.00	36.86	C
ATOM	13152	NZ	LYS	B	327	-37.642	4.436	-22.917	1.00	37.83	N
ATOM	13156	C	LYS	B	327	-42.708	-.752	-21.531	1.00	28.46	C
ATOM	13157	O	LYS	B	327	-42.879	-1.314	-20.461	1.00	28.54	O
ATOM	13159	N	LEU	B	328	-43.696	-.258	-22.268	1.00	27.39	N
ATOM	13160	CA	LEU	B	328	-45.056	-.216	-21.785	1.00	26.79	C
ATOM	13162	CB	LEU	B	328	-45.946	.519	-22.773	1.00	26.68	C
ATOM	13165	CG	LEU	B	328	-47.265	1.034	-22.212	1.00	26.00	C
ATOM	13167	CD1	LEU	B	328	-47.010	1.964	-21.036	1.00	24.88	C
ATOM	13171	CD2	LEU	B	328	-48.051	1.740	-23.308	1.00	24.55	C
ATOM	13175	C	LEU	B	328	-45.560	-1.628	-21.587	1.00	26.70	C
ATOM	13176	O	LEU	B	328	-46.110	-1.983	-20.543	1.00	27.03	O
ATOM	13178	N	CYS	B	329	-45.341	-2.453	-22.591	1.00	26.39	N
ATOM	13179	CA	CYS	B	329	-45.750	-3.843	-22.524	1.00	26.12	C
ATOM	13181	CB	CYS	B	329	-45.382	-4.526	-23.838	1.00	26.47	C

ATOM	13184	SG	CYS B 329	-45.857	-6.236	-23.880	1.00	30.60	S
ATOM	13186	C	CYS B 329	-45.124	-4.568	-21.319	1.00	24.45	C
ATOM	13187	O	CYS B 329	-45.829	-5.105	-20.486	1.00	23.93	O
ATOM	13189	N	PHE B 330	-43.801	-4.553	-21.228	1.00	23.30	N
ATOM	13190	CA	PHE B 330	-43.087	-5.220	-20.141	1.00	22.47	C
ATOM	13192	CB	PHE B 330	-41.575	-5.019	-20.267	1.00	22.42	C
ATOM	13195	CG	PHE B 330	-40.800	-5.518	-19.076	1.00	22.37	C
ATOM	13196	CD1	PHE B 330	-40.502	-6.864	-18.940	1.00	22.47	C
ATOM	13198	CE1	PHE B 330	-39.799	-7.339	-17.836	1.00	21.64	C
ATOM	13200	CZ	PHE B 330	-39.379	-6.472	-16.866	1.00	21.50	C
ATOM	13202	CE2	PHE B 330	-39.659	-5.122	-16.983	1.00	21.92	C
ATOM	13204	CD2	PHE B 330	-40.374	-4.648	-18.084	1.00	22.29	C
ATOM	13206	C	PHE B 330	-43.524	-4.739	-18.765	1.00	21.87	C
ATOM	13207	O	PHE B 330	-43.779	-5.552	-17.875	1.00	21.98	O
ATOM	13209	N	LEU B 331	-43.594	-3.427	-18.573	1.00	20.88	N
ATOM	13210	CA	LEU B 331	-43.997	-2.894	-17.276	1.00	20.25	C
ATOM	13212	CB	LEU B 331	-43.927	-1.354	-17.257	1.00	20.06	C
ATOM	13215	CG	LEU B 331	-44.124	-.628	-15.912	1.00	19.21	C
ATOM	13217	CD1	LEU B 331	-43.542	-1.408	-14.737	1.00	18.78	C
ATOM	13221	CD2	LEU B 331	-43.555	.770	-15.957	1.00	15.75	C
ATOM	13225	C	LEU B 331	-45.398	-3.408	-16.902	1.00	19.99	C
ATOM	13226	O	LEU B 331	-45.606	-3.860	-15.771	1.00	19.67	O
ATOM	13228	N	ALA B 332	-46.336	-3.373	-17.860	1.00	19.56	N
ATOM	13229	CA	ALA B 332	-47.695	-3.882	-17.629	1.00	19.24	C
ATOM	13231	CB	ALA B 332	-48.528	-3.802	-18.891	1.00	18.58	C
ATOM	13235	C	ALA B 332	-47.637	-5.321	-17.108	1.00	19.23	C
ATOM	13236	O	ALA B 332	-48.235	-5.650	-16.079	1.00	19.18	O
ATOM	13238	N	LEU B 333	-46.891	-6.162	-17.816	1.00	19.16	N
ATOM	13239	CA	LEU B 333	-46.727	-7.554	-17.438	1.00	19.08	C
ATOM	13241	CB	LEU B 333	-45.830	-8.260	-18.462	1.00	18.97	C
ATOM	13244	CG	LEU B 333	-45.511	-9.758	-18.307	1.00	19.18	C
ATOM	13246	CD1	LEU B 333	-46.779	-10.611	-18.121	1.00	18.09	C
ATOM	13250	CD2	LEU B 333	-44.676	-10.262	-19.499	1.00	17.17	C
ATOM	13254	C	LEU B 333	-46.111	-7.616	-16.042	1.00	19.16	C
ATOM	13255	O	LEU B 333	-46.629	-8.279	-15.137	1.00	18.97	O
ATOM	13257	N	TYR B 334	-45.013	-6.887	-15.884	1.00	19.23	N
ATOM	13258	CA	TYR B 334	-44.212	-6.906	-14.664	1.00	19.33	C
ATOM	13260	CB	TYR B 334	-43.092	-5.873	-14.785	1.00	19.43	C
ATOM	13263	CG	TYR B 334	-42.190	-5.748	-13.596	1.00	19.44	C
ATOM	13264	CD1	TYR B 334	-41.006	-6.464	-13.524	1.00	20.60	C
ATOM	13266	CE1	TYR B 334	-40.148	-6.347	-12.442	1.00	20.77	C
ATOM	13268	CZ	TYR B 334	-40.474	-5.503	-11.417	1.00	21.51	C
ATOM	13269	OH	TYR B 334	-39.628	-5.387	-10.353	1.00	21.42	O
ATOM	13271	CE2	TYR B 334	-41.646	-4.768	-11.466	1.00	22.15	C
ATOM	13273	CD2	TYR B 334	-42.495	-4.892	-12.565	1.00	20.55	C
ATOM	13275	C	TYR B 334	-45.065	-6.614	-13.443	1.00	19.37	C
ATOM	13276	O	TYR B 334	-44.980	-7.337	-12.441	1.00	19.94	O
ATOM	13278	N	ASN B 335	-45.888	-5.565	-13.527	1.00	18.83	N
ATOM	13279	CA	ASN B 335	-46.784	-5.236	-12.443	1.00	18.37	C
ATOM	13281	CB	ASN B 335	-47.452	-3.896	-12.675	1.00	18.52	C
ATOM	13284	CG	ASN B 335	-46.493	-2.742	-12.518	1.00	19.27	C
ATOM	13285	OD1	ASN B 335	-45.421	-2.906	-11.953	1.00	21.25	O
ATOM	13286	ND2	ASN B 335	-46.872	-1.566	-13.022	1.00	19.24	N
ATOM	13289	C	ASN B 335	-47.812	-6.333	-12.291	1.00	18.22	C
ATOM	13290	O	ASN B 335	-47.966	-6.891	-11.207	1.00	18.53	O
ATOM	13292	N	THR B 336	-48.481	-6.697	-13.379	1.00	17.95	N
ATOM	13293	CA	THR B 336	-49.525	-7.728	-13.301	1.00	17.61	C
ATOM	13295	CB	THR B 336	-49.980	-8.204	-14.676	1.00	17.36	C
ATOM	13297	OG1	THR B 336	-50.249	-7.065	-15.500	1.00	17.02	O
ATOM	13299	CG2	THR B 336	-51.228	-9.044	-14.551	1.00	16.07	C

ATOM	13303	C	THR B 336	-49.065	-8.941	-12.501	1.00	17.75	C
ATOM	13304	O	THR B 336	-49.788	-9.429	-11.621	1.00	18.09	O
ATOM	13306	N	ILE B 337	-47.859	-9.411	-12.785	1.00	17.69	N
ATOM	13307	CA	ILE B 337	-47.384	-10.630	-12.165	1.00	17.74	C
ATOM	13309	CB	ILE B 337	-46.228	-11.245	-12.949	1.00	17.40	C
ATOM	13311	CG1	ILE B 337	-46.795	-11.865	-14.227	1.00	18.44	C
ATOM	13314	CD1	ILE B 337	-45.767	-12.192	-15.300	1.00	19.26	C
ATOM	13318	CG2	ILE B 337	-45.568	-12.320	-12.152	1.00	16.25	C
ATOM	13322	C	ILE B 337	-47.053	-10.395	-10.699	1.00	18.32	C
ATOM	13323	O	ILE B 337	-47.497	-11.175	-9.838	1.00	18.06	O
ATOM	13325	N	ASN B 338	-46.321	-9.311	-10.404	1.00	18.89	N
ATOM	13326	CA	ASN B 338	-45.985	-8.978	-9.007	1.00	19.49	C
ATOM	13328	CB	ASN B 338	-45.189	-7.690	-8.917	1.00	19.43	C
ATOM	13331	CG	ASN B 338	-43.789	-7.836	-9.444	1.00	20.83	C
ATOM	13332	OD1	ASN B 338	-43.292	-8.954	-9.634	1.00	21.52	O
ATOM	13333	ND2	ASN B 338	-43.124	-6.697	-9.676	1.00	22.32	N
ATOM	13336	C	ASN B 338	-47.229	-8.835	-8.146	1.00	20.05	C
ATOM	13337	O	ASN B 338	-47.182	-9.055	-6.953	1.00	19.81	O
ATOM	13339	N	GLU B 339	-48.337	-8.459	-8.770	1.00	20.81	N
ATOM	13340	CA	GLU B 339	-49.589	-8.347	-8.086	1.00	21.81	C
ATOM	13342	CB	GLU B 339	-50.563	-7.544	-8.933	1.00	22.69	C
ATOM	13345	CG	GLU B 339	-51.240	-6.422	-8.148	1.00	27.28	C
ATOM	13348	CD	GLU B 339	-52.571	-5.959	-8.768	1.00	32.79	C
ATOM	13349	OE1	GLU B 339	-52.761	-4.713	-8.897	1.00	36.20	O
ATOM	13350	OE2	GLU B 339	-53.416	-6.837	-9.114	1.00	34.35	O
ATOM	13351	C	GLU B 339	-50.170	-9.732	-7.762	1.00	21.59	C
ATOM	13352	O	GLU B 339	-50.690	-9.946	-6.666	1.00	21.38	O
ATOM	13354	N	ILE B 340	-50.094	-10.676	-8.700	1.00	21.48	N
ATOM	13355	CA	ILE B 340	-50.486	-12.049	-8.384	1.00	21.27	C
ATOM	13357	CB	ILE B 340	-50.437	-12.990	-9.603	1.00	21.17	C
ATOM	13359	CG1	ILE B 340	-51.478	-12.589	-10.643	1.00	20.99	C
ATOM	13362	CD1	ILE B 340	-51.245	-13.202	-12.017	1.00	19.72	C
ATOM	13366	CG2	ILE B 340	-50.702	-14.444	-9.177	1.00	21.09	C
ATOM	13370	C	ILE B 340	-49.557	-12.585	-7.278	1.00	21.40	C
ATOM	13371	O	ILE B 340	-50.023	-13.238	-6.334	1.00	21.74	O
ATOM	13373	N	ALA B 341	-48.256	-12.301	-7.380	1.00	20.90	N
ATOM	13374	CA	ALA B 341	-47.312	-12.721	-6.349	1.00	20.62	C
ATOM	13376	CB	ALA B 341	-45.922	-12.260	-6.694	1.00	20.53	C
ATOM	13380	C	ALA B 341	-47.720	-12.195	-4.969	1.00	20.59	C
ATOM	13381	O	ALA B 341	-47.606	-12.901	-3.950	1.00	20.57	O
ATOM	13383	N	TYR B 342	-48.208	-10.958	-4.949	1.00	20.62	N
ATOM	13384	CA	TYR B 342	-48.600	-10.303	-3.707	1.00	20.59	C
ATOM	13386	CB	TYR B 342	-48.790	-8.792	-3.907	1.00	20.05	C
ATOM	13389	CG	TYR B 342	-49.309	-8.143	-2.674	1.00	18.29	C
ATOM	13390	CD1	TYR B 342	-48.443	-7.683	-1.698	1.00	17.62	C
ATOM	13392	CE1	TYR B 342	-48.914	-7.111	-.518	1.00	17.91	C
ATOM	13394	CZ	TYR B 342	-50.285	-7.005	-.298	1.00	18.18	C
ATOM	13395	OH	TYR B 342	-50.763	-6.433	.877	1.00	16.00	O
ATOM	13397	CE2	TYR B 342	-51.169	-7.467	-1.275	1.00	18.10	C
ATOM	13399	CD2	TYR B 342	-50.670	-8.037	-2.450	1.00	17.67	C
ATOM	13401	C	TYR B 342	-49.863	-10.928	-3.118	1.00	21.48	C
ATOM	13402	O	TYR B 342	-49.960	-11.089	-1.913	1.00	20.91	O
ATOM	13404	N	ASP B 343	-50.832	-11.273	-3.959	1.00	22.93	N
ATOM	13405	CA	ASP B 343	-52.048	-11.925	-3.470	1.00	24.30	C
ATOM	13407	CB	ASP B 343	-53.010	-12.282	-4.608	1.00	24.59	C
ATOM	13410	CG	ASP B 343	-53.577	-11.060	-5.334	1.00	26.66	C
ATOM	13411	OD1	ASP B 343	-53.619	-9.936	-4.756	1.00	28.24	O
ATOM	13412	OD2	ASP B 343	-54.004	-11.247	-6.506	1.00	29.56	O
ATOM	13413	C	ASP B 343	-51.660	-13.208	-2.761	1.00	24.97	C
ATOM	13414	O	ASP B 343	-52.128	-13.495	-1.658	1.00	24.95	O

ATOM	13416	N	ASN B 344	-50.792	-13.973	-3.413	1.00	25.90	N
ATOM	13417	CA	ASN B 344	-50.322	-15.252	-2.884	1.00	26.59	C
ATOM	13419	CB	ASN B 344	-49.623	-16.043	-3.978	1.00	26.74	C
ATOM	13422	CG	ASN B 344	-50.594	-16.639	-4.933	1.00	28.21	C
ATOM	13423	OD1	ASN B 344	-51.080	-17.737	-4.688	1.00	32.10	O
ATOM	13424	ND2	ASN B 344	-50.923	-15.916	-6.016	1.00	28.40	N
ATOM	13427	C	ASN B 344	-49.421	-15.146	-1.655	1.00	26.78	C
ATOM	13428	O	ASN B 344	-49.424	-16.053	-.821	1.00	27.06	O
ATOM	13430	N	LEU B 345	-48.647	-14.067	-1.533	1.00	26.74	N
ATOM	13431	CA	LEU B 345	-47.970	-13.793	-.261	1.00	26.59	C
ATOM	13433	CB	LEU B 345	-47.005	-12.618	-.396	1.00	26.25	C
ATOM	13436	CG	LEU B 345	-46.046	-12.376	.764	1.00	24.65	C
ATOM	13438	CD1	LEU B 345	-45.258	-13.621	1.115	1.00	22.50	C
ATOM	13442	CD2	LEU B 345	-45.119	-11.251	.387	1.00	23.54	C
ATOM	13446	C	LEU B 345	-49.004	-13.503	.840	1.00	27.00	C
ATOM	13447	O	LEU B 345	-48.903	-14.014	1.947	1.00	26.93	O
ATOM	13449	N	LYS B 346	-50.008	-12.697	.518	1.00	27.54	N
ATOM	13450	CA	LYS B 346	-51.012	-12.300	1.491	1.00	27.96	C
ATOM	13452	CB	LYS B 346	-51.998	-11.288	.889	1.00	28.13	C
ATOM	13455	CG	LYS B 346	-52.822	-10.527	1.926	1.00	28.70	C
ATOM	13458	CD	LYS B 346	-53.781	-9.526	1.281	1.00	29.73	C
ATOM	13461	CE	LYS B 346	-55.205	-10.063	1.133	1.00	31.17	C
ATOM	13464	NZ	LYS B 346	-55.839	-9.661	-.171	1.00	32.83	N
ATOM	13468	C	LYS B 346	-51.788	-13.498	1.984	1.00	28.31	C
ATOM	13469	O	LYS B 346	-52.074	-13.600	3.177	1.00	28.52	O
ATOM	13471	N	ASP B 347	-52.149	-14.396	1.075	1.00	28.56	N
ATOM	13472	CA	ASP B 347	-53.102	-15.441	1.427	1.00	29.11	C
ATOM	13474	CB	ASP B 347	-54.068	-15.736	.268	1.00	29.50	C
ATOM	13477	CG	ASP B 347	-54.902	-14.500	-.150	1.00	31.04	C
ATOM	13478	OD1	ASP B 347	-54.961	-13.500	.616	1.00	31.35	O
ATOM	13479	OD2	ASP B 347	-55.496	-14.534	-1.261	1.00	33.39	O
ATOM	13480	C	ASP B 347	-52.403	-16.699	1.891	1.00	28.79	C
ATOM	13481	O	ASP B 347	-52.824	-17.321	2.854	1.00	29.02	O
ATOM	13483	N	LYS B 348	-51.326	-17.068	1.222	1.00	28.74	N
ATOM	13484	CA	LYS B 348	-50.612	-18.293	1.564	1.00	28.84	C
ATOM	13486	CB	LYS B 348	-50.290	-19.101	.299	1.00	29.24	C
ATOM	13489	CG	LYS B 348	-51.513	-19.470	-.576	1.00	30.97	C
ATOM	13492	CD	LYS B 348	-51.051	-20.300	-1.803	1.00	33.66	C
ATOM	13495	CE	LYS B 348	-52.033	-20.244	-2.989	1.00	34.52	C
ATOM	13498	NZ	LYS B 348	-53.445	-20.542	-2.607	1.00	35.61	N
ATOM	13502	C	LYS B 348	-49.330	-18.029	2.358	1.00	27.97	C
ATOM	13503	O	LYS B 348	-48.763	-18.948	2.925	1.00	27.88	O
ATOM	13505	N	GLY B 349	-48.874	-16.785	2.403	1.00	27.21	N
ATOM	13506	CA	GLY B 349	-47.625	-16.473	3.075	1.00	26.87	C
ATOM	13509	C	GLY B 349	-46.437	-17.227	2.514	1.00	26.71	C
ATOM	13510	O	GLY B 349	-45.576	-17.669	3.271	1.00	26.65	O
ATOM	13512	N	GLU B 350	-46.409	-17.390	1.191	1.00	26.61	N
ATOM	13513	CA	GLU B 350	-45.263	-17.953	.470	1.00	26.43	C
ATOM	13515	CB	GLU B 350	-45.624	-19.310	-.128	1.00	26.85	C
ATOM	13518	CG	GLU B 350	-45.685	-20.453	.888	1.00	29.63	C
ATOM	13521	CD	GLU B 350	-44.304	-21.036	1.227	1.00	34.21	C
ATOM	13522	OE1	GLU B 350	-43.417	-21.045	.328	1.00	36.93	O
ATOM	13523	OE2	GLU B 350	-44.106	-21.498	2.384	1.00	35.89	O
ATOM	13524	C	GLU B 350	-44.900	-16.981	-.637	1.00	25.49	C
ATOM	13525	O	GLU B 350	-45.774	-16.312	-1.181	1.00	25.35	O
ATOM	13527	N	ASN B 351	-43.616	-16.866	-.955	1.00	24.74	N
ATOM	13528	CA	ASN B 351	-43.206	-16.046	-2.090	1.00	24.47	C
ATOM	13530	CB	ASN B 351	-41.851	-15.369	-1.867	1.00	24.94	C
ATOM	13533	CG	ASN B 351	-41.428	-14.489	-3.062	1.00	26.98	C
ATOM	13534	OD1	ASN B 351	-41.976	-14.605	-4.170	1.00	28.96	O

ATOM	13535	ND2 ASN B 351	-40.457	-13.601	-2.834	1.00	29.55	N
ATOM	13538	C ASN B 351	-43.095	-16.918	-3.303	1.00	23.33	C
ATOM	13539	O ASN B 351	-42.120	-17.641	-3.429	1.00	23.47	O
ATOM	13541	N ILE B 352	-44.060	-16.823	-4.208	1.00	22.09	N
ATOM	13542	CA ILE B 352	-44.048	-17.644	-5.410	1.00	21.22	C
ATOM	13544	CB ILE B 352	-45.452	-18.213	-5.675	1.00	21.23	C
ATOM	13546	CG1 ILE B 352	-46.464	-17.080	-5.913	1.00	21.11	C
ATOM	13549	CD1 ILE B 352	-47.346	-17.305	-7.126	1.00	20.33	C
ATOM	13553	CG2 ILE B 352	-45.875	-19.116	-4.519	1.00	19.79	C
ATOM	13557	C ILE B 352	-43.508	-16.947	-6.689	1.00	20.88	C
ATOM	13558	O ILE B 352	-43.486	-17.553	-7.758	1.00	20.76	O
ATOM	13560	N LEU B 353	-43.030	-15.705	-6.574	1.00	20.44	N
ATOM	13561	CA LEU B 353	-42.603	-14.907	-7.750	1.00	19.84	C
ATOM	13563	CB LEU B 353	-42.055	-13.527	-7.325	1.00	19.67	C
ATOM	13566	CG LEU B 353	-42.003	-12.339	-8.311	1.00	18.47	C
ATOM	13568	CD1 LEU B 353	-43.306	-12.083	-9.008	1.00	17.03	C
ATOM	13572	CD2 LEU B 353	-41.604	-11.064	-7.597	1.00	17.08	C
ATOM	13576	C LEU B 353	-41.587	-15.628	-8.633	1.00	19.81	C
ATOM	13577	O LEU B 353	-41.720	-15.596	-9.851	1.00	19.55	O
ATOM	13579	N PRO B 354	-40.578	-16.299	-8.024	1.00	20.01	N
ATOM	13580	CA PRO B 354	-39.571	-17.010	-8.822	1.00	19.73	C
ATOM	13582	CB PRO B 354	-38.781	-17.796	-7.779	1.00	19.55	C
ATOM	13585	CG PRO B 354	-38.928	-17.037	-6.542	1.00	19.61	C
ATOM	13588	CD PRO B 354	-40.289	-16.413	-6.578	1.00	19.90	C
ATOM	13591	C PRO B 354	-40.187	-17.971	-9.804	1.00	19.87	C
ATOM	13592	O PRO B 354	-39.693	-18.103	-10.918	1.00	20.08	O
ATOM	13593	N TYR B 355	-41.264	-18.636	-9.387	1.00	20.14	N
ATOM	13594	CA TYR B 355	-41.909	-19.665	-10.216	1.00	20.24	C
ATOM	13596	CB TYR B 355	-42.878	-20.520	-9.394	1.00	20.40	C
ATOM	13599	CG TYR B 355	-42.189	-21.153	-8.214	1.00	21.82	C
ATOM	13600	CD1 TYR B 355	-41.119	-21.999	-8.412	1.00	22.86	C
ATOM	13602	CE1 TYR B 355	-40.457	-22.564	-7.353	1.00	25.07	C
ATOM	13604	CZ TYR B 355	-40.850	-22.285	-6.053	1.00	25.58	C
ATOM	13605	OH TYR B 355	-40.153	-22.887	-5.020	1.00	27.33	O
ATOM	13607	CE2 TYR B 355	-41.918	-21.433	-5.815	1.00	23.93	C
ATOM	13609	CD2 TYR B 355	-42.577	-20.867	-6.900	1.00	23.10	C
ATOM	13611	C TYR B 355	-42.625	-19.001	-11.361	1.00	19.67	C
ATOM	13612	O TYR B 355	-42.455	-19.405	-12.510	1.00	19.75	O
ATOM	13614	N LEU B 356	-43.385	-17.954	-11.030	1.00	18.97	N
ATOM	13615	CA LEU B 356	-44.180	-17.205	-12.004	1.00	18.26	C
ATOM	13617	CB LEU B 356	-45.039	-16.157	-11.293	1.00	18.03	C
ATOM	13620	CG LEU B 356	-46.056	-16.727	-10.312	1.00	18.13	C
ATOM	13622	CD1 LEU B 356	-46.691	-15.613	-9.500	1.00	19.37	C
ATOM	13626	CD2 LEU B 356	-47.106	-17.521	-11.046	1.00	17.93	C
ATOM	13630	C LEU B 356	-43.298	-16.536	-13.056	1.00	17.76	C
ATOM	13631	O LEU B 356	-43.598	-16.567	-14.259	1.00	16.93	O
ATOM	13633	N THR B 357	-42.201	-15.949	-12.600	1.00	17.56	N
ATOM	13634	CA THR B 357	-41.335	-15.215	-13.501	1.00	17.63	C
ATOM	13636	CB THR B 357	-40.438	-14.236	-12.759	1.00	17.69	C
ATOM	13638	OG1 THR B 357	-39.530	-14.974	-11.934	1.00	18.39	O
ATOM	13640	CG2 THR B 357	-41.288	-13.242	-11.927	1.00	16.47	C
ATOM	13644	C THR B 357	-40.477	-16.145	-14.356	1.00	17.55	C
ATOM	13645	O THR B 357	-40.154	-15.793	-15.498	1.00	17.17	O
ATOM	13647	N LYS B 358	-40.127	-17.321	-13.824	1.00	17.37	N
ATOM	13648	CA LYS B 358	-39.462	-18.337	-14.649	1.00	17.55	C
ATOM	13650	CB LYS B 358	-39.030	-19.542	-13.827	1.00	17.87	C
ATOM	13653	CG LYS B 358	-38.450	-20.718	-14.640	1.00	18.75	C
ATOM	13656	CD LYS B 358	-37.013	-20.497	-15.112	1.00	20.18	C
ATOM	13659	CE LYS B 358	-36.361	-21.843	-15.457	1.00	21.49	C
ATOM	13662	NZ LYS B 358	-35.189	-21.702	-16.361	1.00	22.89	N

ATOM	13666	C	LYS B 358	-40.403	-18.777	-15.754	1.00	17.43	C
ATOM	13667	O	LYS B 358	-40.031	-18.775	-16.926	1.00	17.52	O
ATOM	13669	N	ALA B 359	-41.633	-19.120	-15.384	1.00	17.42	N
ATOM	13670	CA	ALA B 359	-42.654	-19.471	-16.359	1.00	17.52	C
ATOM	13672	CB	ALA B 359	-44.011	-19.434	-15.738	1.00	17.16	C
ATOM	13676	C	ALA B 359	-42.585	-18.521	-17.531	1.00	18.18	C
ATOM	13677	O	ALA B 359	-42.513	-18.956	-18.675	1.00	18.50	O
ATOM	13679	N	TRP B 360	-42.556	-17.223	-17.245	1.00	19.14	N
ATOM	13680	CA	TRP B 360	-42.526	-16.206	-18.301	1.00	19.77	C
ATOM	13682	CB	TRP B 360	-42.922	-14.837	-17.746	1.00	20.15	C
ATOM	13685	CG	TRP B 360	-44.377	-14.653	-17.813	1.00	20.50	C
ATOM	13686	CD1	TRP B 360	-45.259	-14.729	-16.788	1.00	22.15	C
ATOM	13688	NE1	TRP B 360	-46.536	-14.533	-17.248	1.00	22.84	N
ATOM	13690	CE2	TRP B 360	-46.487	-14.348	-18.603	1.00	21.68	C
ATOM	13691	CD2	TRP B 360	-45.134	-14.418	-18.987	1.00	20.85	C
ATOM	13692	CE3	TRP B 360	-44.802	-14.253	-20.334	1.00	21.21	C
ATOM	13694	CZ3	TRP B 360	-45.821	-14.025	-21.245	1.00	20.75	C
ATOM	13696	CH2	TRP B 360	-47.166	-13.961	-20.827	1.00	21.21	C
ATOM	13698	CZ2	TRP B 360	-47.515	-14.123	-19.516	1.00	20.97	C
ATOM	13700	C	TRP B 360	-41.215	-16.110	-19.081	1.00	19.96	C
ATOM	13701	O	TRP B 360	-41.242	-15.791	-20.271	1.00	19.94	O
ATOM	13703	N	ALA B 361	-40.089	-16.368	-18.420	1.00	20.14	N
ATOM	13704	CA	ALA B 361	-38.790	-16.437	-19.102	1.00	20.41	C
ATOM	13706	CB	ALA B 361	-37.676	-16.581	-18.102	1.00	20.67	C
ATOM	13710	C	ALA B 361	-38.763	-17.618	-20.041	1.00	20.55	C
ATOM	13711	O	ALA B 361	-38.330	-17.510	-21.187	1.00	20.76	O
ATOM	13713	N	ASP B 362	-39.230	-18.749	-19.534	1.00	20.57	N
ATOM	13714	CA	ASP B 362	-39.360	-19.965	-20.326	1.00	20.97	C
ATOM	13716	CB	ASP B 362	-39.984	-21.082	-19.470	1.00	21.39	C
ATOM	13719	CG	ASP B 362	-38.938	-21.902	-18.693	1.00	22.80	C
ATOM	13720	OD1	ASP B 362	-38.836	-23.109	-18.970	1.00	28.57	O
ATOM	13721	OD2	ASP B 362	-38.210	-21.384	-17.825	1.00	22.78	O
ATOM	13722	C	ASP B 362	-40.179	-19.748	-21.618	1.00	20.66	C
ATOM	13723	O	ASP B 362	-39.740	-20.141	-22.701	1.00	20.69	O
ATOM	13725	N	LEU B 363	-41.346	-19.110	-21.501	1.00	20.22	N
ATOM	13726	CA	LEU B 363	-42.245	-18.885	-22.648	1.00	19.82	C
ATOM	13728	CB	LEU B 363	-43.592	-18.327	-22.164	1.00	19.11	C
ATOM	13731	CG	LEU B 363	-44.644	-17.945	-23.209	1.00	17.02	C
ATOM	13733	CD1	LEU B 363	-44.829	-19.033	-24.214	1.00	16.20	C
ATOM	13737	CD2	LEU B 363	-45.979	-17.618	-22.586	1.00	14.21	C
ATOM	13741	C	LEU B 363	-41.633	-17.950	-23.706	1.00	20.76	C
ATOM	13742	O	LEU B 363	-41.656	-18.235	-24.910	1.00	21.00	O
ATOM	13744	N	CYS B 364	-41.088	-16.826	-23.253	1.00	21.42	N
ATOM	13745	CA	CYS B 364	-40.427	-15.884	-24.147	1.00	21.67	C
ATOM	13747	CB	CYS B 364	-39.969	-14.639	-23.376	1.00	21.86	C
ATOM	13750	SG	CYS B 364	-41.327	-13.638	-22.722	1.00	23.36	S
ATOM	13752	C	CYS B 364	-39.229	-16.547	-24.835	1.00	21.52	C
ATOM	13753	O	CYS B 364	-39.010	-16.323	-26.037	1.00	21.67	O
ATOM	13755	N	ASN B 365	-38.454	-17.352	-24.095	1.00	20.83	N
ATOM	13756	CA	ASN B 365	-37.306	-18.025	-24.713	1.00	20.65	C
ATOM	13758	CB	ASN B 365	-36.425	-18.756	-23.691	1.00	20.74	C
ATOM	13761	CG	ASN B 365	-35.330	-17.865	-23.100	1.00	20.53	C
ATOM	13762	OD1	ASN B 365	-34.553	-17.229	-23.830	1.00	18.26	O
ATOM	13763	ND2	ASN B 365	-35.252	-17.841	-21.763	1.00	20.69	N
ATOM	13766	C	ASN B 365	-37.756	-18.987	-25.811	1.00	20.43	C
ATOM	13767	O	ASN B 365	-36.998	-19.253	-26.759	1.00	20.11	O
ATOM	13769	N	ALA B 366	-38.987	-19.492	-25.667	1.00	20.01	N
ATOM	13770	CA	ALA B 366	-39.618	-20.340	-26.665	1.00	19.70	C
ATOM	13772	CB	ALA B 366	-40.766	-21.113	-26.045	1.00	19.53	C
ATOM	13776	C	ALA B 366	-40.099	-19.494	-27.838	1.00	19.74	C

ATOM 13777	O	ALA B 366	-39.846	-19.824	-28.989	1.00	19.61	O
ATOM 13779	N	PHE B 367	-40.791	-18.399	-27.553	1.00	19.99	N
ATOM 13780	CA	PHE B 367	-41.138	-17.446	-28.606	1.00	20.36	C
ATOM 13782	CB	PHE B 367	-41.805	-16.194	-28.026	1.00	20.62	C
ATOM 13785	CG	PHE B 367	-43.200	-16.403	-27.537	1.00	21.75	C
ATOM 13786	CD1	PHE B 367	-44.136	-17.071	-28.315	1.00	23.00	C
ATOM 13788	CE1	PHE B 367	-45.435	-17.242	-27.874	1.00	23.60	C
ATOM 13790	CZ	PHE B 367	-45.819	-16.731	-26.648	1.00	23.64	C
ATOM 13792	CE2	PHE B 367	-44.901	-16.049	-25.871	1.00	23.56	C
ATOM 13794	CD2	PHE B 367	-43.599	-15.883	-26.316	1.00	22.91	C
ATOM 13796	C	PHE B 367	-39.906	-16.986	-29.393	1.00	20.19	C
ATOM 13797	O	PHE B 367	-39.934	-16.903	-30.615	1.00	19.98	O
ATOM 13799	N	LEU B 368	-38.839	-16.659	-28.675	1.00	20.24	N
ATOM 13800	CA	LEU B 368	-37.629	-16.161	-29.298	1.00	20.31	C
ATOM 13802	CB	LEU B 368	-36.569	-15.900	-28.232	1.00	20.27	C
ATOM 13805	CG	LEU B 368	-35.235	-15.313	-28.694	1.00	20.12	C
ATOM 13807	CD1	LEU B 368	-35.417	-14.069	-29.592	1.00	20.05	C
ATOM 13811	CD2	LEU B 368	-34.409	-14.987	-27.465	1.00	18.60	C
ATOM 13815	C	LEU B 368	-37.122	-17.179	-30.295	1.00	20.53	C
ATOM 13816	O	LEU B 368	-36.844	-16.844	-31.446	1.00	20.49	O
ATOM 13818	N	GLN B 369	-37.032	-18.426	-29.832	1.00	20.80	N
ATOM 13819	CA	GLN B 369	-36.565	-19.557	-30.638	1.00	20.90	C
ATOM 13821	CB	GLN B 369	-36.614	-20.857	-29.815	1.00	20.98	C
ATOM 13824	CG	GLN B 369	-36.343	-22.159	-30.576	1.00	20.75	C
ATOM 13827	CD	GLN B 369	-34.946	-22.238	-31.126	1.00	20.29	C
ATOM 13828	OE1	GLN B 369	-34.718	-21.984	-32.305	1.00	20.17	O
ATOM 13829	NE2	GLN B 369	-33.996	-22.581	-30.272	1.00	20.33	N
ATOM 13832	C	GLN B 369	-37.365	-19.718	-31.910	1.00	21.04	C
ATOM 13833	O	GLN B 369	-36.803	-20.097	-32.929	1.00	20.82	O
ATOM 13835	N	GLU B 370	-38.668	-19.453	-31.863	1.00	21.80	N
ATOM 13836	CA	GLU B 370	-39.480	-19.512	-33.087	1.00	22.73	C
ATOM 13838	CB	GLU B 370	-40.996	-19.518	-32.801	1.00	23.18	C
ATOM 13841	CG	GLU B 370	-41.449	-20.608	-31.809	1.00	25.92	C
ATOM 13844	CD	GLU B 370	-42.927	-21.064	-31.957	1.00	29.74	C
ATOM 13845	OE1	GLU B 370	-43.840	-20.191	-32.101	1.00	30.13	O
ATOM 13846	OE2	GLU B 370	-43.156	-22.316	-31.892	1.00	31.94	O
ATOM 13847	C	GLU B 370	-39.079	-18.358	-34.016	1.00	22.62	C
ATOM 13848	O	GLU B 370	-38.933	-18.565	-35.221	1.00	22.53	O
ATOM 13850	N	ALA B 371	-38.864	-17.166	-33.456	1.00	22.65	N
ATOM 13851	CA	ALA B 371	-38.427	-16.023	-34.254	1.00	22.95	C
ATOM 13853	CB	ALA B 371	-38.375	-14.753	-33.424	1.00	22.97	C
ATOM 13857	C	ALA B 371	-37.070	-16.296	-34.891	1.00	23.13	C
ATOM 13858	O	ALA B 371	-36.900	-16.063	-36.091	1.00	23.38	O
ATOM 13860	N	LYS B 372	-36.113	-16.804	-34.107	1.00	23.15	N
ATOM 13861	CA	LYS B 372	-34.776	-17.102	-34.640	1.00	23.10	C
ATOM 13863	CB	LYS B 372	-33.808	-17.633	-33.575	1.00	22.85	C
ATOM 13866	CG	LYS B 372	-33.414	-16.616	-32.514	1.00	24.28	C
ATOM 13869	CD	LYS B 372	-31.973	-16.806	-31.969	1.00	26.34	C
ATOM 13872	CE	LYS B 372	-31.919	-17.407	-30.544	1.00	27.84	C
ATOM 13875	NZ	LYS B 372	-31.673	-16.405	-29.458	1.00	27.93	N
ATOM 13879	C	LYS B 372	-34.887	-18.093	-35.790	1.00	23.16	C
ATOM 13880	O	LYS B 372	-34.231	-17.920	-36.804	1.00	23.73	O
ATOM 13882	N	TRP B 373	-35.720	-19.122	-35.658	1.00	23.11	N
ATOM 13883	CA	TRP B 373	-35.844	-20.094	-36.742	1.00	22.99	C
ATOM 13885	CB	TRP B 373	-36.746	-21.274	-36.370	1.00	22.71	C
ATOM 13888	CG	TRP B 373	-36.081	-22.343	-35.562	1.00	20.71	C
ATOM 13889	CD1	TRP B 373	-34.770	-22.740	-35.623	1.00	19.89	C
ATOM 13891	NE1	TRP B 373	-34.540	-23.757	-34.734	1.00	19.05	N
ATOM 13893	CE2	TRP B 373	-35.715	-24.050	-34.090	1.00	18.91	C
ATOM 13894	CD2	TRP B 373	-36.707	-23.184	-34.599	1.00	18.41	C

ATOM	13895	CE3 TRP B 373	-38.007	-23.288	-34.105	1.00	15.66	C
ATOM	13897	CZ3 TRP B 373	-38.271	-24.216	-33.130	1.00	15.41	C
ATOM	13899	CH2 TRP B 373	-37.276	-25.069	-32.647	1.00	16.33	C
ATOM	13901	CZ2 TRP B 373	-35.994	-25.005	-33.112	1.00	17.53	C
ATOM	13903	C TRP B 373	-36.374	-19.416	-37.991	1.00	23.59	C
ATOM	13904	O TRP B 373	-35.878	-19.653	-39.096	1.00	23.60	O
ATOM	13906	N LEU B 374	-37.374	-18.559	-37.801	1.00	24.32	N
ATOM	13907	CA LEU B 374	-38.008	-17.859	-38.915	1.00	24.66	C
ATOM	13909	CB LEU B 374	-39.181	-17.019	-38.428	1.00	24.85	C
ATOM	13912	CG LEU B 374	-40.112	-16.449	-39.496	1.00	25.45	C
ATOM	13914	CD1 LEU B 374	-41.002	-17.535	-40.070	1.00	24.83	C
ATOM	13918	CD2 LEU B 374	-40.959	-15.293	-38.885	1.00	26.92	C
ATOM	13922	C LEU B 374	-37.009	-16.975	-39.618	1.00	24.64	C
ATOM	13923	O LEU B 374	-36.934	-16.992	-40.830	1.00	24.33	O
ATOM	13925	N TYR B 375	-36.232	-16.220	-38.850	1.00	25.09	N
ATOM	13926	CA TYR B 375	-35.283	-15.274	-39.433	1.00	25.71	C
ATOM	13928	CB TYR B 375	-34.534	-14.479	-38.349	1.00	25.91	C
ATOM	13931	CG TYR B 375	-33.536	-13.468	-38.892	1.00	27.11	C
ATOM	13932	CD1 TYR B 375	-33.950	-12.200	-39.305	1.00	27.98	C
ATOM	13934	CE1 TYR B 375	-33.039	-11.271	-39.811	1.00	28.61	C
ATOM	13936	CZ TYR B 375	-31.693	-11.605	-39.904	1.00	29.13	C
ATOM	13937	OH TYR B 375	-30.790	-10.685	-40.397	1.00	29.83	O
ATOM	13939	CE2 TYR B 375	-31.256	-12.863	-39.501	1.00	28.79	C
ATOM	13941	CD2 TYR B 375	-32.177	-13.783	-38.995	1.00	28.03	C
ATOM	13943	C TYR B 375	-34.305	-16.036	-40.301	1.00	25.64	C
ATOM	13944	O TYR B 375	-34.201	-15.788	-41.503	1.00	25.74	O
ATOM	13946	N ASN B 376	-33.635	-17.000	-39.685	1.00	25.73	N
ATOM	13947	CA ASN B 376	-32.599	-17.777	-40.346	1.00	25.78	C
ATOM	13949	CB ASN B 376	-31.739	-18.497	-39.303	1.00	25.73	C
ATOM	13952	CG ASN B 376	-31.155	-17.565	-38.273	1.00	25.13	C
ATOM	13953	OD1 ASN B 376	-30.491	-16.590	-38.601	1.00	25.53	O
ATOM	13954	ND2 ASN B 376	-31.389	-17.874	-37.013	1.00	24.64	N
ATOM	13957	C ASN B 376	-33.126	-18.824	-41.332	1.00	26.04	C
ATOM	13958	O ASN B 376	-32.343	-19.644	-41.807	1.00	26.36	O
ATOM	13960	N LYS B 377	-34.428	-18.824	-41.632	1.00	26.16	N
ATOM	13961	CA LYS B 377	-35.012	-19.811	-42.552	1.00	26.27	C
ATOM	13963	CB LYS B 377	-34.575	-19.538	-44.010	1.00	26.51	C
ATOM	13966	CG LYS B 377	-35.470	-18.561	-44.802	1.00	28.01	C
ATOM	13969	CD LYS B 377	-34.629	-17.592	-45.667	1.00	29.99	C
ATOM	13972	CE LYS B 377	-35.493	-16.605	-46.465	1.00	31.02	C
ATOM	13975	NZ LYS B 377	-36.420	-15.819	-45.589	1.00	32.17	N
ATOM	13979	C LYS B 377	-34.647	-21.243	-42.150	1.00	25.93	C
ATOM	13980	O LYS B 377	-34.380	-22.084	-43.003	1.00	25.82	O
ATOM	13982	N SER B 378	-34.637	-21.522	-40.852	1.00	25.65	N
ATOM	13983	CA SER B 378	-34.324	-22.860	-40.385	1.00	25.54	C
ATOM	13985	CB SER B 378	-34.132	-22.877	-38.882	1.00	25.41	C
ATOM	13988	OG SER B 378	-32.905	-22.269	-38.570	1.00	25.67	O
ATOM	13990	C SER B 378	-35.412	-23.838	-40.790	1.00	25.62	C
ATOM	13991	O SER B 378	-36.495	-23.435	-41.199	1.00	25.70	O
ATOM	13993	N THR B 379	-35.108	-25.129	-40.690	1.00	25.66	N
ATOM	13994	CA THR B 379	-36.038	-26.175	-41.097	1.00	25.55	C
ATOM	13996	CB THR B 379	-35.788	-26.635	-42.546	1.00	25.59	C
ATOM	13998	OG1 THR B 379	-34.391	-26.899	-42.736	1.00	25.46	O
ATOM	14000	CG2 THR B 379	-36.255	-25.574	-43.528	1.00	25.88	C
ATOM	14004	C THR B 379	-35.886	-27.348	-40.160	1.00	25.45	C
ATOM	14005	O THR B 379	-35.372	-28.393	-40.548	1.00	25.70	O
ATOM	14007	N PRO B 380	-36.335	-27.179	-38.915	1.00	25.37	N
ATOM	14008	CA PRO B 380	-36.144	-28.182	-37.876	1.00	25.22	C
ATOM	14010	CB PRO B 380	-36.500	-27.437	-36.583	1.00	25.19	C
ATOM	14013	CG PRO B 380	-36.643	-26.021	-36.955	1.00	25.83	C

ATOM	14016	CD PRO B 380	-37.025	-25.998	-38.390	1.00	25.71	C
ATOM	14019	C PRO B 380	-37.052	-29.379	-38.009	1.00	24.82	C
ATOM	14020	O PRO B 380	-38.143	-29.280	-38.575	1.00	24.84	O
ATOM	14021	N THR B 381	-36.607	-30.493	-37.436	1.00	24.42	N
ATOM	14022	CA THR B 381	-37.386	-31.712	-37.429	1.00	23.99	C
ATOM	14024	CB THR B 381	-36.614	-32.873	-36.789	1.00	23.99	C
ATOM	14026	OG1 THR B 381	-36.339	-32.575	-35.418	1.00	24.02	O
ATOM	14028	CG2 THR B 381	-35.302	-33.117	-37.524	1.00	23.76	C
ATOM	14032	C THR B 381	-38.649	-31.461	-36.636	1.00	23.78	C
ATOM	14033	O THR B 381	-38.700	-30.548	-35.813	1.00	23.89	O
ATOM	14035	N PHE B 382	-39.672	-32.267	-36.889	1.00	23.53	N
ATOM	14036	CA PHE B 382	-40.905	-32.204	-36.105	1.00	22.97	C
ATOM	14038	CB PHE B 382	-41.870	-33.305	-36.527	1.00	22.89	C
ATOM	14041	CG PHE B 382	-43.079	-33.394	-35.655	1.00	22.26	C
ATOM	14042	CD1 PHE B 382	-44.189	-32.601	-35.910	1.00	21.65	C
ATOM	14044	CE1 PHE B 382	-45.300	-32.672	-35.099	1.00	20.87	C
ATOM	14046	CZ PHE B 382	-45.310	-33.534	-34.016	1.00	20.36	C
ATOM	14048	CE2 PHE B 382	-44.204	-34.317	-33.744	1.00	20.52	C
ATOM	14050	CD2 PHE B 382	-43.098	-34.243	-34.556	1.00	21.17	C
ATOM	14052	C PHE B 382	-40.677	-32.329	-34.605	1.00	22.69	C
ATOM	14053	O PHE B 382	-41.341	-31.657	-33.834	1.00	22.44	O
ATOM	14055	N ASP B 383	-39.764	-33.205	-34.191	1.00	22.52	N
ATOM	14056	CA ASP B 383	-39.491	-33.375	-32.765	1.00	22.50	C
ATOM	14058	CB ASP B 383	-38.531	-34.525	-32.527	1.00	22.43	C
ATOM	14061	CG ASP B 383	-39.174	-35.863	-32.711	1.00	22.71	C
ATOM	14062	OD1 ASP B 383	-40.359	-35.945	-33.074	1.00	22.89	O
ATOM	14063	OD2 ASP B 383	-38.470	-36.856	-32.492	1.00	25.17	O
ATOM	14064	C ASP B 383	-38.954	-32.116	-32.091	1.00	22.43	C
ATOM	14065	O ASP B 383	-39.438	-31.758	-31.014	1.00	22.40	O
ATOM	14067	N ASP B 384	-37.966	-31.464	-32.713	1.00	22.37	N
ATOM	14068	CA ASP B 384	-37.422	-30.191	-32.213	1.00	22.61	C
ATOM	14070	CB ASP B 384	-36.317	-29.648	-33.124	1.00	22.82	C
ATOM	14073	CG ASP B 384	-34.963	-30.208	-32.807	1.00	23.40	C
ATOM	14074	OD1 ASP B 384	-34.847	-30.973	-31.832	1.00	25.33	O
ATOM	14075	OD2 ASP B 384	-34.009	-29.891	-33.545	1.00	24.26	O
ATOM	14076	C ASP B 384	-38.482	-29.113	-32.130	1.00	22.56	C
ATOM	14077	O ASP B 384	-38.598	-28.429	-31.108	1.00	22.48	O
ATOM	14079	N TYR B 385	-39.227	-28.963	-33.227	1.00	22.36	N
ATOM	14080	CA TYR B 385	-40.170	-27.870	-33.402	1.00	22.18	C
ATOM	14082	CB TYR B 385	-40.738	-27.867	-34.812	1.00	22.05	C
ATOM	14085	CG TYR B 385	-41.818	-26.834	-35.019	1.00	22.60	C
ATOM	14086	CD1 TYR B 385	-41.493	-25.553	-35.437	1.00	23.64	C
ATOM	14088	CE1 TYR B 385	-42.462	-24.580	-35.630	1.00	25.00	C
ATOM	14090	CZ TYR B 385	-43.788	-24.882	-35.410	1.00	26.86	C
ATOM	14091	OH TYR B 385	-44.734	-23.883	-35.622	1.00	28.48	O
ATOM	14093	CE2 TYR B 385	-44.149	-26.171	-34.980	1.00	25.95	C
ATOM	14095	CD2 TYR B 385	-43.158	-27.133	-34.790	1.00	23.94	C
ATOM	14097	C TYR B 385	-41.319	-27.986	-32.446	1.00	22.13	C
ATOM	14098	O TYR B 385	-41.703	-27.017	-31.808	1.00	22.95	O
ATOM	14100	N PHE B 386	-41.899	-29.169	-32.380	1.00	21.80	N
ATOM	14101	CA PHE B 386	-43.048	-29.394	-31.525	1.00	21.62	C
ATOM	14103	CB PHE B 386	-43.624	-30.795	-31.768	1.00	21.73	C
ATOM	14106	CG PHE B 386	-44.834	-31.104	-30.952	1.00	21.27	C
ATOM	14107	CD1 PHE B 386	-46.034	-30.447	-31.196	1.00	22.55	C
ATOM	14109	CE1 PHE B 386	-47.164	-30.732	-30.448	1.00	23.14	C
ATOM	14111	CZ PHE B 386	-47.095	-31.699	-29.446	1.00	23.20	C
ATOM	14113	CE2 PHE B 386	-45.896	-32.358	-29.208	1.00	21.57	C
ATOM	14115	CD2 PHE B 386	-44.781	-32.057	-29.959	1.00	20.29	C
ATOM	14117	C PHE B 386	-42.630	-29.226	-30.073	1.00	21.34	C
ATOM	14118	O PHE B 386	-43.353	-28.639	-29.284	1.00	21.44	O

ATOM	14120	N	GLY B 387	-41.447	-29.725	-29.736	1.00	20.84	N
ATOM	14121	CA	GLY B 387	-40.935	-29.618	-28.387	1.00	20.57	C
ATOM	14124	C	GLY B 387	-40.854	-28.186	-27.916	1.00	20.26	C
ATOM	14125	O	GLY B 387	-40.930	-27.908	-26.724	1.00	20.86	O
ATOM	14127	N	ASN B 388	-40.691	-27.268	-28.852	1.00	19.81	N
ATOM	14128	CA	ASN B 388	-40.758	-25.853	-28.536	1.00	19.51	C
ATOM	14130	CB	ASN B 388	-39.877	-25.086	-29.508	1.00	19.42	C
ATOM	14133	CG	ASN B 388	-39.593	-23.700	-29.045	1.00	19.26	C
ATOM	14134	OD1	ASN B 388	-38.916	-23.513	-28.038	1.00	18.93	O
ATOM	14135	ND2	ASN B 388	-40.105	-22.709	-29.774	1.00	17.99	N
ATOM	14138	C	ASN B 388	-42.200	-25.320	-28.604	1.00	19.41	C
ATOM	14139	O	ASN B 388	-42.591	-24.447	-27.833	1.00	19.58	O
ATOM	14141	N	ALA B 389	-42.989	-25.848	-29.532	1.00	19.11	N
ATOM	14142	CA	ALA B 389	-44.313	-25.307	-29.808	1.00	18.87	C
ATOM	14144	CB	ALA B 389	-44.897	-25.949	-31.060	1.00	18.87	C
ATOM	14148	C	ALA B 389	-45.264	-25.460	-28.634	1.00	18.51	C
ATOM	14149	O	ALA B 389	-45.916	-24.511	-28.262	1.00	18.33	O
ATOM	14151	N	TRP B 390	-45.344	-26.652	-28.054	1.00	18.73	N
ATOM	14152	CA	TRP B 390	-46.258	-26.887	-26.935	1.00	18.74	C
ATOM	14154	CB	TRP B 390	-46.389	-28.377	-26.569	1.00	18.75	C
ATOM	14157	CG	TRP B 390	-45.166	-29.122	-26.019	1.00	18.46	C
ATOM	14158	CD1	TRP B 390	-44.390	-30.001	-26.705	1.00	18.75	C
ATOM	14160	NE1	TRP B 390	-43.417	-30.519	-25.897	1.00	17.56	N
ATOM	14162	CE2	TRP B 390	-43.559	-30.005	-24.642	1.00	16.67	C
ATOM	14163	CD2	TRP B 390	-44.661	-29.124	-24.674	1.00	17.50	C
ATOM	14164	CE3	TRP B 390	-45.008	-28.448	-23.505	1.00	18.38	C
ATOM	14166	CZ3	TRP B 390	-44.251	-28.685	-22.347	1.00	18.34	C
ATOM	14168	CH2	TRP B 390	-43.164	-29.573	-22.361	1.00	16.75	C
ATOM	14170	CZ2	TRP B 390	-42.805	-30.234	-23.495	1.00	15.75	C
ATOM	14172	C	TRP B 390	-45.861	-26.066	-25.727	1.00	18.91	C
ATOM	14173	O	TRP B 390	-46.707	-25.691	-24.919	1.00	19.02	O
ATOM	14175	N	LYS B 391	-44.570	-25.787	-25.612	1.00	18.95	N
ATOM	14176	CA	LYS B 391	-44.085	-24.810	-24.643	1.00	18.85	C
ATOM	14178	CB	LYS B 391	-42.544	-24.888	-24.508	1.00	19.52	C
ATOM	14181	CG	LYS B 391	-42.023	-25.373	-23.143	1.00	21.17	C
ATOM	14184	CD	LYS B 391	-40.516	-25.689	-23.191	1.00	22.87	C
ATOM	14187	CE	LYS B 391	-40.251	-27.159	-23.508	1.00	23.91	C
ATOM	14190	NZ	LYS B 391	-38.938	-27.354	-24.191	1.00	25.44	N
ATOM	14194	C	LYS B 391	-44.537	-23.400	-25.032	1.00	17.57	C
ATOM	14195	O	LYS B 391	-44.937	-22.629	-24.180	1.00	17.66	O
ATOM	14197	N	SER B 392	-44.490	-23.068	-26.314	1.00	16.65	N
ATOM	14198	CA	SER B 392	-44.879	-21.718	-26.764	1.00	16.27	C
ATOM	14200	CB	SER B 392	-44.241	-21.375	-28.119	1.00	16.22	C
ATOM	14203	OG	SER B 392	-44.937	-21.969	-29.207	1.00	15.94	O
ATOM	14205	C	SER B 392	-46.384	-21.504	-26.866	1.00	15.87	C
ATOM	14206	O	SER B 392	-46.825	-20.395	-27.093	1.00	15.56	O
ATOM	14208	N	SER B 393	-47.167	-22.568	-26.724	1.00	15.92	N
ATOM	14209	CA	SER B 393	-48.629	-22.474	-26.718	1.00	15.80	C
ATOM	14211	CB	SER B 393	-49.240	-23.867	-26.630	1.00	15.78	C
ATOM	14214	OG	SER B 393	-49.025	-24.426	-25.348	1.00	15.00	O
ATOM	14216	C	SER B 393	-49.097	-21.646	-25.533	1.00	15.94	C
ATOM	14217	O	SER B 393	-50.115	-20.948	-25.599	1.00	15.60	O
ATOM	14219	N	SER B 394	-48.296	-21.740	-24.471	1.00	16.14	N
ATOM	14220	CA	SER B 394	-48.487	-21.091	-23.177	1.00	16.31	C
ATOM	14222	CB	SER B 394	-49.062	-19.662	-23.271	1.00	16.18	C
ATOM	14225	OG	SER B 394	-50.472	-19.649	-23.358	1.00	16.57	O
ATOM	14227	C	SER B 394	-49.316	-22.006	-22.297	1.00	16.34	C
ATOM	14228	O	SER B 394	-49.822	-21.595	-21.261	1.00	16.47	O
ATOM	14230	N	GLY B 395	-49.403	-23.268	-22.702	1.00	16.55	N
ATOM	14231	CA	GLY B 395	-50.103	-24.282	-21.927	1.00	16.79	C

ATOM	14234	C	GLY B 395	-49.543	-24.391	-20.534	1.00	16.80	C
ATOM	14235	O	GLY B 395	-50.222	-24.083	-19.556	1.00	16.74	O
ATOM	14237	N	PRO B 396	-48.291	-24.824	-20.430	1.00	17.07	N
ATOM	14238	CA	PRO B 396	-47.698	-24.930	-19.103	1.00	17.18	C
ATOM	14240	CB	PRO B 396	-46.273	-25.399	-19.386	1.00	17.09	C
ATOM	14243	CG	PRO B 396	-46.087	-25.237	-20.874	1.00	17.77	C
ATOM	14246	CD	PRO B 396	-47.424	-25.390	-21.470	1.00	17.20	C
ATOM	14249	C	PRO B 396	-47.707	-23.615	-18.313	1.00	17.17	C
ATOM	14250	O	PRO B 396	-47.921	-23.644	-17.089	1.00	17.10	O
ATOM	14251	N	LEU B 397	-47.499	-22.475	-18.983	1.00	17.04	N
ATOM	14252	CA	LEU B 397	-47.513	-21.183	-18.261	1.00	16.77	C
ATOM	14254	CB	LEU B 397	-47.116	-19.969	-19.135	1.00	16.84	C
ATOM	14257	CG	LEU B 397	-47.145	-18.576	-18.458	1.00	16.91	C
ATOM	14259	CD1	LEU B 397	-46.577	-18.641	-17.096	1.00	18.36	C
ATOM	14263	CD2	LEU B 397	-46.373	-17.515	-19.196	1.00	16.88	C
ATOM	14267	C	LEU B 397	-48.894	-20.988	-17.682	1.00	16.24	C
ATOM	14268	O	LEU B 397	-49.051	-20.710	-16.494	1.00	16.09	O
ATOM	14270	N	GLN B 398	-49.902	-21.173	-18.515	1.00	15.72	N
ATOM	14271	CA	GLN B 398	-51.262	-21.093	-18.024	1.00	15.49	C
ATOM	14273	CB	GLN B 398	-52.267	-21.387	-19.120	1.00	15.51	C
ATOM	14276	CG	GLN B 398	-52.371	-20.275	-20.118	1.00	16.12	C
ATOM	14279	CD	GLN B 398	-53.436	-20.545	-21.114	1.00	17.86	C
ATOM	14280	OE1	GLN B 398	-54.509	-21.028	-20.757	1.00	20.99	O
ATOM	14281	NE2	GLN B 398	-53.170	-20.239	-22.374	1.00	18.31	N
ATOM	14284	C	GLN B 398	-51.471	-22.040	-16.873	1.00	15.03	C
ATOM	14285	O	GLN B 398	-51.974	-21.638	-15.843	1.00	15.21	O
ATOM	14287	N	LEU B 399	-51.065	-23.291	-17.022	1.00	14.72	N
ATOM	14288	CA	LEU B 399	-51.361	-24.254	-15.978	1.00	14.65	C
ATOM	14290	CB	LEU B 399	-51.201	-25.688	-16.475	1.00	14.44	C
ATOM	14293	CG	LEU B 399	-52.250	-26.191	-17.478	1.00	14.28	C
ATOM	14295	CD1	LEU B 399	-51.907	-27.633	-17.846	1.00	15.50	C
ATOM	14299	CD2	LEU B 399	-53.713	-26.077	-16.986	1.00	10.96	C
ATOM	14303	C	LEU B 399	-50.554	-23.995	-14.704	1.00	14.88	C
ATOM	14304	O	LEU B 399	-51.100	-24.161	-13.618	1.00	15.21	O
ATOM	14306	N	ILE B 400	-49.291	-23.562	-14.810	1.00	14.91	N
ATOM	14307	CA	ILE B 400	-48.532	-23.174	-13.607	1.00	14.89	C
ATOM	14309	CB	ILE B 400	-47.158	-22.574	-13.907	1.00	15.19	C
ATOM	14311	CG1	ILE B 400	-46.189	-23.674	-14.353	1.00	16.76	C
ATOM	14314	CD1	ILE B 400	-44.777	-23.162	-14.716	1.00	18.19	C
ATOM	14318	CG2	ILE B 400	-46.603	-21.906	-12.665	1.00	13.65	C
ATOM	14322	C	ILE B 400	-49.288	-22.137	-12.819	1.00	14.77	C
ATOM	14323	O	ILE B 400	-49.485	-22.302	-11.632	1.00	15.36	O
ATOM	14325	N	PHE B 401	-49.717	-21.071	-13.486	1.00	14.59	N
ATOM	14326	CA	PHE B 401	-50.491	-20.001	-12.844	1.00	14.22	C
ATOM	14328	CB	PHE B 401	-50.825	-18.900	-13.845	1.00	14.12	C
ATOM	14331	CG	PHE B 401	-49.803	-17.790	-13.872	1.00	13.83	C
ATOM	14332	CD1	PHE B 401	-50.012	-16.623	-13.173	1.00	12.77	C
ATOM	14334	CE1	PHE B 401	-49.074	-15.629	-13.189	1.00	13.09	C
ATOM	14336	CZ	PHE B 401	-47.906	-15.781	-13.902	1.00	12.96	C
ATOM	14338	CE2	PHE B 401	-47.684	-16.928	-14.592	1.00	12.86	C
ATOM	14340	CD2	PHE B 401	-48.622	-17.932	-14.574	1.00	13.30	C
ATOM	14342	C	PHE B 401	-51.765	-20.478	-12.212	1.00	14.26	C
ATOM	14343	O	PHE B 401	-52.184	-19.944	-11.207	1.00	13.90	O
ATOM	14345	N	ALA B 402	-52.377	-21.480	-12.828	1.00	15.09	N
ATOM	14346	CA	ALA B 402	-53.641	-22.054	-12.364	1.00	15.76	C
ATOM	14348	CB	ALA B 402	-54.264	-22.899	-13.454	1.00	15.57	C
ATOM	14352	C	ALA B 402	-53.417	-22.896	-11.129	1.00	16.63	C
ATOM	14353	O	ALA B 402	-54.284	-22.970	-10.259	1.00	17.00	O
ATOM	14355	N	TYR B 403	-52.253	-23.544	-11.060	1.00	17.50	N
ATOM	14356	CA	TYR B 403	-51.885	-24.317	-9.894	1.00	17.86	C

ATOM	14358	CB	TYR	B	403	-50.486	-24.911	-10.038	1.00	17.73	C
ATOM	14361	CG	TYR	B	403	-50.006	-25.576	-8.764	1.00	18.24	C
ATOM	14362	CD1	TYR	B	403	-50.401	-26.867	-8.436	1.00	18.12	C
ATOM	14364	CE1	TYR	B	403	-49.971	-27.476	-7.273	1.00	18.07	C
ATOM	14366	CZ	TYR	B	403	-49.145	-26.790	-6.409	1.00	18.62	C
ATOM	14367	OH	TYR	B	403	-48.727	-27.383	-5.244	1.00	17.70	O
ATOM	14369	CE2	TYR	B	403	-48.748	-25.499	-6.703	1.00	18.93	C
ATOM	14371	CD2	TYR	B	403	-49.177	-24.901	-7.876	1.00	18.88	C
ATOM	14373	C	TYR	B	403	-51.966	-23.461	-8.630	1.00	18.48	C
ATOM	14374	O	TYR	B	403	-52.494	-23.908	-7.616	1.00	18.66	O
ATOM	14376	N	PHE	B	404	-51.468	-22.232	-8.682	1.00	18.91	N
ATOM	14377	CA	PHE	B	404	-51.400	-21.430	-7.468	1.00	19.53	C
ATOM	14379	CB	PHE	B	404	-50.395	-20.325	-7.644	1.00	19.35	C
ATOM	14382	CG	PHE	B	404	-49.014	-20.808	-7.799	1.00	18.79	C
ATOM	14383	CD1	PHE	B	404	-48.311	-21.246	-6.701	1.00	17.67	C
ATOM	14385	CE1	PHE	B	404	-47.001	-21.683	-6.832	1.00	18.44	C
ATOM	14387	CZ	PHE	B	404	-46.384	-21.688	-8.078	1.00	18.44	C
ATOM	14389	CE2	PHE	B	404	-47.083	-21.245	-9.190	1.00	18.79	C
ATOM	14391	CD2	PHE	B	404	-48.396	-20.805	-9.047	1.00	18.90	C
ATOM	14393	C	PHE	B	404	-52.733	-20.817	-7.072	1.00	20.49	C
ATOM	14394	O	PHE	B	404	-52.925	-20.389	-5.924	1.00	19.95	O
ATOM	14396	N	ALA	B	405	-53.636	-20.756	-8.043	1.00	21.92	N
ATOM	14397	CA	ALA	B	405	-54.918	-20.106	-7.872	1.00	23.09	C
ATOM	14399	CB	ALA	B	405	-55.333	-19.426	-9.167	1.00	23.15	C
ATOM	14403	C	ALA	B	405	-55.959	-21.113	-7.446	1.00	24.23	C
ATOM	14404	O	ALA	B	405	-57.003	-20.726	-6.925	1.00	24.33	O
ATOM	14406	N	VAL	B	406	-55.657	-22.399	-7.662	1.00	25.80	N
ATOM	14407	CA	VAL	B	406	-56.569	-23.511	-7.355	1.00	26.88	C
ATOM	14409	CB	VAL	B	406	-56.640	-24.514	-8.512	1.00	26.64	C
ATOM	14411	CG1	VAL	B	406	-57.132	-25.854	-8.012	1.00	27.24	C
ATOM	14415	CG2	VAL	B	406	-57.547	-23.986	-9.593	1.00	26.45	C
ATOM	14419	C	VAL	B	406	-56.161	-24.266	-6.094	1.00	27.91	C
ATOM	14420	O	VAL	B	406	-56.932	-24.345	-5.155	1.00	28.28	O
ATOM	14422	N	VAL	B	407	-54.954	-24.825	-6.087	1.00	29.25	N
ATOM	14423	CA	VAL	B	407	-54.443	-25.573	-4.941	1.00	30.16	C
ATOM	14425	CB	VAL	B	407	-53.128	-26.276	-5.279	1.00	30.04	C
ATOM	14427	CG1	VAL	B	407	-52.482	-26.831	-4.032	1.00	30.42	C
ATOM	14431	CG2	VAL	B	407	-53.378	-27.376	-6.280	1.00	30.14	C
ATOM	14435	C	VAL	B	407	-54.208	-24.646	-3.755	1.00	31.34	C
ATOM	14436	O	VAL	B	407	-53.535	-23.618	-3.876	1.00	31.46	O
ATOM	14438	N	GLN	B	408	-54.753	-25.032	-2.604	1.00	32.63	N
ATOM	14439	CA	GLN	B	408	-54.727	-24.191	-1.417	1.00	33.51	C
ATOM	14441	CB	GLN	B	408	-55.891	-24.572	-.514	1.00	33.84	C
ATOM	14444	CG	GLN	B	408	-56.161	-23.548	.577	1.00	35.49	C
ATOM	14447	CD	GLN	B	408	-57.623	-23.161	.656	1.00	37.61	C
ATOM	14448	OE1	GLN	B	408	-58.519	-24.006	.501	1.00	38.23	O
ATOM	14449	NE2	GLN	B	408	-57.876	-21.874	.893	1.00	38.64	N
ATOM	14452	C	GLN	B	408	-53.399	-24.274	-.647	1.00	33.72	C
ATOM	14453	O	GLN	B	408	-52.852	-23.253	-.204	1.00	33.52	O
ATOM	14455	N	ASN	B	409	-52.889	-25.491	-.490	1.00	33.95	N
ATOM	14456	CA	ASN	B	409	-51.642	-25.707	.226	1.00	34.16	C
ATOM	14458	CB	ASN	B	409	-51.865	-26.665	1.391	1.00	34.25	C
ATOM	14461	CG	ASN	B	409	-52.756	-26.069	2.459	1.00	34.39	C
ATOM	14462	OD1	ASN	B	409	-52.269	-25.577	3.480	1.00	34.40	O
ATOM	14463	ND2	ASN	B	409	-54.070	-26.089	2.222	1.00	34.12	N
ATOM	14466	C	ASN	B	409	-50.582	-26.244	-.709	1.00	34.12	C
ATOM	14467	O	ASN	B	409	-50.578	-27.422	-1.046	1.00	34.41	O
ATOM	14469	N	ILE	B	410	-49.681	-25.369	-1.127	1.00	34.07	N
ATOM	14470	CA	ILE	B	410	-48.699	-25.718	-2.138	1.00	34.08	C
ATOM	14472	CB	ILE	B	410	-48.138	-24.455	-2.840	1.00	34.15	C

ATOM	14474	CG1 ILE B 410	-47.274	-23.610	-1.891	1.00	34.26	C
ATOM	14477	CD1 ILE B 410	-47.216	-22.139	-2.249	1.00	34.10	C
ATOM	14481	CG2 ILE B 410	-49.290	-23.634	-3.404	1.00	34.49	C
ATOM	14485	C ILE B 410	-47.586	-26.553	-1.533	1.00	33.93	C
ATOM	14486	O ILE B 410	-47.181	-26.317	-.405	1.00	33.80	O
ATOM	14488	N LYS B 411	-47.123	-27.546	-2.285	1.00	34.12	N
ATOM	14489	CA LYS B 411	-46.012	-28.395	-1.874	1.00	34.45	C
ATOM	14491	CB LYS B 411	-46.414	-29.873	-1.907	1.00	34.67	C
ATOM	14494	CG LYS B 411	-47.850	-30.130	-1.460	1.00	35.79	C
ATOM	14497	CD LYS B 411	-48.102	-31.586	-1.052	1.00	37.60	C
ATOM	14500	CE LYS B 411	-49.450	-31.728	-.309	1.00	38.80	C
ATOM	14503	NZ LYS B 411	-49.568	-32.992	.488	1.00	39.20	N
ATOM	14507	C LYS B 411	-44.843	-28.132	-2.810	1.00	34.33	C
ATOM	14508	O LYS B 411	-45.038	-27.956	-4.006	1.00	34.14	O
ATOM	14510	N LYS B 412	-43.631	-28.102	-2.265	1.00	34.50	N
ATOM	14511	CA LYS B 412	-42.460	-27.688	-3.042	1.00	34.69	C
ATOM	14513	CB LYS B 412	-41.242	-27.398	-2.154	1.00	35.13	C
ATOM	14516	CG LYS B 412	-41.205	-25.960	-1.613	1.00	36.90	C
ATOM	14519	CD LYS B 412	-40.079	-25.749	-.588	1.00	38.57	C
ATOM	14522	CE LYS B 412	-40.546	-24.852	.555	1.00	39.55	C
ATOM	14525	NZ LYS B 412	-39.491	-24.644	1.587	1.00	40.83	N
ATOM	14529	C LYS B 412	-42.075	-28.687	-4.103	1.00	34.13	C
ATOM	14530	O LYS B 412	-41.468	-28.308	-5.093	1.00	34.21	O
ATOM	14532	N GLU B 413	-42.408	-29.958	-3.910	1.00	33.54	N
ATOM	14533	CA GLU B 413	-42.095	-30.949	-4.933	1.00	33.19	C
ATOM	14535	CB GLU B 413	-41.886	-32.335	-4.330	1.00	33.44	C
ATOM	14538	CG GLU B 413	-43.127	-33.037	-3.807	1.00	34.37	C
ATOM	14541	CD GLU B 413	-42.834	-34.490	-3.499	1.00	35.56	C
ATOM	14542	OE1 GLU B 413	-42.662	-35.264	-4.471	1.00	35.30	O
ATOM	14543	OE2 GLU B 413	-42.751	-34.847	-2.297	1.00	36.69	O
ATOM	14544	C GLU B 413	-43.157	-30.965	-6.029	1.00	32.47	C
ATOM	14545	O GLU B 413	-42.846	-31.232	-7.193	1.00	32.28	O
ATOM	14547	N GLU B 414	-44.403	-30.676	-5.652	1.00	31.58	N
ATOM	14548	CA GLU B 414	-45.482	-30.512	-6.614	1.00	30.99	C
ATOM	14550	CB GLU B 414	-46.781	-30.101	-5.927	1.00	30.98	C
ATOM	14553	CG GLU B 414	-47.732	-31.245	-5.642	1.00	31.82	C
ATOM	14556	CD GLU B 414	-49.100	-30.774	-5.138	1.00	34.29	C
ATOM	14557	OE1 GLU B 414	-49.258	-29.581	-4.777	1.00	35.53	O
ATOM	14558	OE2 GLU B 414	-50.036	-31.603	-5.101	1.00	36.44	O
ATOM	14559	C GLU B 414	-45.104	-29.455	-7.628	1.00	30.59	C
ATOM	14560	O GLU B 414	-45.169	-29.687	-8.828	1.00	30.37	O
ATOM	14562	N ILE B 415	-44.684	-28.295	-7.140	1.00	30.49	N
ATOM	14563	CA ILE B 415	-44.367	-27.177	-8.028	1.00	30.40	C
ATOM	14565	CB ILE B 415	-44.412	-25.797	-7.320	1.00	30.40	C
ATOM	14567	CG1 ILE B 415	-43.235	-25.589	-6.388	1.00	30.36	C
ATOM	14570	CD1 ILE B 415	-43.373	-24.313	-5.611	1.00	31.01	C
ATOM	14574	CG2 ILE B 415	-45.686	-25.644	-6.523	1.00	30.85	C
ATOM	14578	C ILE B 415	-43.037	-27.347	-8.730	1.00	30.12	C
ATOM	14579	O ILE B 415	-42.870	-26.865	-9.840	1.00	30.42	O
ATOM	14581	N GLU B 416	-42.095	-28.030	-8.099	1.00	29.75	N
ATOM	14582	CA GLU B 416	-40.799	-28.249	-8.719	1.00	29.70	C
ATOM	14584	CB GLU B 416	-39.825	-28.777	-7.690	1.00	30.02	C
ATOM	14587	CG GLU B 416	-38.386	-28.426	-7.948	1.00	31.20	C
ATOM	14590	CD GLU B 416	-37.523	-28.839	-6.776	1.00	32.91	C
ATOM	14591	OE1 GLU B 416	-38.008	-28.723	-5.632	1.00	32.37	O
ATOM	14592	OE2 GLU B 416	-36.377	-29.294	-6.994	1.00	35.31	O
ATOM	14593	C GLU B 416	-40.917	-29.222	-9.890	1.00	29.24	C
ATOM	14594	O GLU B 416	-40.121	-29.177	-10.835	1.00	28.86	O
ATOM	14596	N ASN B 417	-41.915	-30.097	-9.819	1.00	28.81	N
ATOM	14597	CA ASN B 417	-42.252	-30.960	-10.941	1.00	28.66	C

ATOM	14599	CB	ASN B 417	-43.165	-32.105	-10.503	1.00	28.72	C
ATOM	14602	CG	ASN B 417	-42.379	-33.285	-9.973	1.00	29.78	C
ATOM	14603	OD1	ASN B 417	-41.887	-34.108	-10.744	1.00	30.39	O
ATOM	14604	ND2	ASN B 417	-42.223	-33.358	-8.653	1.00	31.45	N
ATOM	14607	C	ASN B 417	-42.888	-30.183	-12.067	1.00	28.32	C
ATOM	14608	O	ASN B 417	-42.611	-30.456	-13.232	1.00	27.95	O
ATOM	14610	N	LEU B 418	-43.740	-29.217	-11.713	1.00	28.23	N
ATOM	14611	CA	LEU B 418	-44.349	-28.313	-12.697	1.00	27.92	C
ATOM	14613	CB	LEU B 418	-45.298	-27.320	-12.023	1.00	27.48	C
ATOM	14616	CG	LEU B 418	-46.636	-27.896	-11.553	1.00	27.02	C
ATOM	14618	CD1	LEU B 418	-47.393	-26.882	-10.691	1.00	26.70	C
ATOM	14622	CD2	LEU B 418	-47.499	-28.355	-12.720	1.00	25.90	C
ATOM	14626	C	LEU B 418	-43.275	-27.567	-13.491	1.00	28.20	C
ATOM	14627	O	LEU B 418	-43.310	-27.541	-14.733	1.00	27.90	O
ATOM	14629	N	GLN B 419	-42.308	-26.995	-12.770	1.00	28.43	N
ATOM	14630	CA	GLN B 419	-41.170	-26.315	-13.395	1.00	28.72	C
ATOM	14632	CB	GLN B 419	-40.223	-25.746	-12.347	1.00	28.81	C
ATOM	14635	CG	GLN B 419	-40.592	-24.332	-11.946	1.00	30.14	C
ATOM	14638	CD	GLN B 419	-39.535	-23.671	-11.092	1.00	32.05	C
ATOM	14639	OE1	GLN B 419	-39.246	-22.477	-11.257	1.00	33.88	O
ATOM	14640	NE2	GLN B 419	-38.948	-24.437	-10.170	1.00	31.84	N
ATOM	14643	C	GLN B 419	-40.390	-27.179	-14.370	1.00	28.83	C
ATOM	14644	O	GLN B 419	-39.922	-26.668	-15.386	1.00	28.72	O
ATOM	14646	N	LYS B 420	-40.265	-28.475	-14.067	1.00	29.21	N
ATOM	14647	CA	LYS B 420	-39.613	-29.451	-14.968	1.00	29.26	C
ATOM	14649	CB	LYS B 420	-38.924	-30.555	-14.143	1.00	29.42	C
ATOM	14652	CG	LYS B 420	-37.800	-30.052	-13.207	1.00	30.58	C
ATOM	14655	CD	LYS B 420	-37.373	-31.113	-12.151	1.00	32.35	C
ATOM	14658	CE	LYS B 420	-36.572	-30.506	-10.955	1.00	33.05	C
ATOM	14661	NZ	LYS B 420	-36.520	-31.368	-9.702	1.00	32.62	N
ATOM	14665	C	LYS B 420	-40.579	-30.057	-16.021	1.00	28.96	C
ATOM	14666	O	LYS B 420	-40.216	-30.965	-16.753	1.00	28.66	O
ATOM	14668	N	TYR B 421	-41.804	-29.543	-16.088	1.00	29.01	N
ATOM	14669	CA	TYR B 421	-42.777	-29.893	-17.134	1.00	29.16	C
ATOM	14671	CB	TYR B 421	-42.209	-29.602	-18.534	1.00	29.48	C
ATOM	14674	CG	TYR B 421	-41.968	-28.127	-18.773	1.00	31.04	C
ATOM	14675	CD1	TYR B 421	-43.018	-27.219	-18.716	1.00	32.54	C
ATOM	14677	CE1	TYR B 421	-42.815	-25.878	-18.914	1.00	33.35	C
ATOM	14679	CZ	TYR B 421	-41.556	-25.413	-19.193	1.00	34.56	C
ATOM	14680	OH	TYR B 421	-41.377	-24.066	-19.384	1.00	37.99	O
ATOM	14682	CE2	TYR B 421	-40.490	-26.278	-19.267	1.00	33.56	C
ATOM	14684	CD2	TYR B 421	-40.702	-27.638	-19.051	1.00	32.75	C
ATOM	14686	C	TYR B 421	-43.345	-31.311	-17.042	1.00	28.60	C
ATOM	14687	O	TYR B 421	-43.395	-32.046	-18.025	1.00	28.49	O
ATOM	14689	N	HIS B 422	-43.808	-31.662	-15.846	1.00	28.18	N
ATOM	14690	CA	HIS B 422	-44.507	-32.919	-15.594	1.00	27.66	C
ATOM	14692	CB	HIS B 422	-45.082	-32.914	-14.180	1.00	27.73	C
ATOM	14695	CG	HIS B 422	-45.500	-34.264	-13.685	1.00	28.14	C
ATOM	14696	ND1	HIS B 422	-44.597	-35.195	-13.217	1.00	29.15	N
ATOM	14698	CE1	HIS B 422	-45.247	-36.277	-12.829	1.00	29.50	C
ATOM	14700	NE2	HIS B 422	-46.539	-36.078	-13.022	1.00	28.69	N
ATOM	14702	CD2	HIS B 422	-46.724	-34.827	-13.552	1.00	27.24	C
ATOM	14704	C	HIS B 422	-45.650	-33.164	-16.575	1.00	27.29	C
ATOM	14705	O	HIS B 422	-46.463	-32.269	-16.853	1.00	26.86	O
ATOM	14707	N	ASP B 423	-45.712	-34.402	-17.060	1.00	26.96	N
ATOM	14708	CA	ASP B 423	-46.772	-34.883	-17.952	1.00	26.78	C
ATOM	14710	CB	ASP B 423	-46.792	-36.413	-17.939	1.00	27.01	C
ATOM	14713	CG	ASP B 423	-45.594	-37.018	-18.646	1.00	28.21	C
ATOM	14714	OD1	ASP B 423	-45.100	-36.378	-19.602	1.00	31.33	O
ATOM	14715	OD2	ASP B 423	-45.153	-38.130	-18.262	1.00	28.39	O

ATOM	14716	C	ASP B 423	-48.187	-34.383	-17.651	1.00	26.21	C
ATOM	14717	O	ASP B 423	-49.015	-34.336	-18.556	1.00	26.75	O
ATOM	14719	N	ILE B 424	-48.463	-34.042	-16.391	1.00	25.11	N
ATOM	14720	CA	ILE B 424	-49.773	-33.553	-15.956	1.00	23.82	C
ATOM	14722	CB	ILE B 424	-49.798	-33.310	-14.439	1.00	23.67	C
ATOM	14724	CG1	ILE B 424	-51.196	-32.990	-13.938	1.00	22.97	C
ATOM	14727	CD1	ILE B 424	-51.220	-32.693	-12.465	1.00	22.06	C
ATOM	14731	CG2	ILE B 424	-48.859	-32.171	-14.065	1.00	24.35	C
ATOM	14735	C	ILE B 424	-50.102	-32.265	-16.669	1.00	23.08	C
ATOM	14736	O	ILE B 424	-51.265	-32.025	-17.023	1.00	22.99	O
ATOM	14738	N	ILE B 425	-49.082	-31.436	-16.886	1.00	22.22	N
ATOM	14739	CA	ILE B 425	-49.282	-30.199	-17.633	1.00	21.61	C
ATOM	14741	CB	ILE B 425	-48.796	-28.959	-16.848	1.00	21.18	C
ATOM	14743	CG1	ILE B 425	-47.279	-28.853	-16.804	1.00	19.47	C
ATOM	14746	CD1	ILE B 425	-46.850	-27.581	-16.140	1.00	18.00	C
ATOM	14750	CG2	ILE B 425	-49.357	-28.977	-15.428	1.00	20.75	C
ATOM	14754	C	ILE B 425	-48.665	-30.255	-19.033	1.00	21.73	C
ATOM	14755	O	ILE B 425	-49.100	-29.523	-19.912	1.00	21.75	O
ATOM	14757	N	SER B 426	-47.686	-31.134	-19.258	1.00	21.74	N
ATOM	14758	CA	SER B 426	-47.016	-31.174	-20.560	1.00	21.53	C
ATOM	14760	CB	SER B 426	-45.639	-31.876	-20.516	1.00	21.66	C
ATOM	14763	OG	SER B 426	-45.724	-33.281	-20.405	1.00	22.19	O
ATOM	14765	C	SER B 426	-47.922	-31.798	-21.582	1.00	21.19	C
ATOM	14766	O	SER B 426	-47.950	-31.358	-22.721	1.00	21.37	O
ATOM	14768	N	ARG B 427	-48.692	-32.800	-21.177	1.00	20.91	N
ATOM	14769	CA	ARG B 427	-49.532	-33.504	-22.143	1.00	20.78	C
ATOM	14771	CB	ARG B 427	-50.038	-34.840	-21.616	1.00	20.94	C
ATOM	14774	CG	ARG B 427	-49.006	-35.887	-21.874	1.00	22.76	C
ATOM	14777	CD	ARG B 427	-49.158	-37.101	-21.025	1.00	26.75	C
ATOM	14780	NE	ARG B 427	-47.960	-37.927	-21.185	1.00	29.36	N
ATOM	14782	CZ	ARG B 427	-47.574	-38.887	-20.352	1.00	30.25	C
ATOM	14783	NH1	ARG B 427	-48.293	-39.179	-19.266	1.00	30.27	N
ATOM	14786	NH2	ARG B 427	-46.447	-39.547	-20.611	1.00	31.25	N
ATOM	14789	C	ARG B 427	-50.646	-32.650	-22.668	1.00	20.07	C
ATOM	14790	O	ARG B 427	-50.724	-32.466	-23.866	1.00	20.13	O
ATOM	14792	N	PRO B 428	-51.482	-32.091	-21.786	1.00	19.38	N
ATOM	14793	CA	PRO B 428	-52.533	-31.202	-22.264	1.00	19.05	C
ATOM	14795	CB	PRO B 428	-53.046	-30.547	-20.990	1.00	19.13	C
ATOM	14798	CG	PRO B 428	-52.696	-31.476	-19.918	1.00	19.55	C
ATOM	14801	CD	PRO B 428	-51.443	-32.161	-20.319	1.00	19.20	C
ATOM	14804	C	PRO B 428	-52.016	-30.135	-23.227	1.00	18.84	C
ATOM	14805	O	PRO B 428	-52.688	-29.802	-24.205	1.00	18.73	O
ATOM	14806	N	SER B 429	-50.828	-29.608	-22.964	1.00	18.57	N
ATOM	14807	CA	SER B 429	-50.208	-28.679	-23.899	1.00	18.64	C
ATOM	14809	CB	SER B 429	-48.960	-28.086	-23.291	1.00	18.46	C
ATOM	14812	OG	SER B 429	-49.141	-27.980	-21.911	1.00	19.08	O
ATOM	14814	C	SER B 429	-49.877	-29.315	-25.260	1.00	18.52	C
ATOM	14815	O	SER B 429	-49.886	-28.633	-26.276	1.00	19.02	O
ATOM	14817	N	HIS B 430	-49.579	-30.606	-25.297	1.00	18.29	N
ATOM	14818	CA	HIS B 430	-49.457	-31.273	-26.583	1.00	18.25	C
ATOM	14820	CB	HIS B 430	-49.085	-32.765	-26.467	1.00	18.40	C
ATOM	14823	CG	HIS B 430	-47.753	-33.024	-25.806	1.00	19.41	C
ATOM	14824	ND1	HIS B 430	-46.776	-32.055	-25.657	1.00	19.55	N
ATOM	14826	CE1	HIS B 430	-45.729	-32.576	-25.043	1.00	17.77	C
ATOM	14828	NE2	HIS B 430	-45.981	-33.849	-24.800	1.00	18.52	N
ATOM	14830	CD2	HIS B 430	-47.232	-34.160	-25.277	1.00	19.19	C
ATOM	14832	C	HIS B 430	-50.789	-31.105	-27.306	1.00	17.93	C
ATOM	14833	O	HIS B 430	-50.816	-30.599	-28.417	1.00	18.54	O
ATOM	14835	N	ILE B 431	-51.897	-31.479	-26.665	1.00	17.34	N
ATOM	14836	CA	ILE B 431	-53.217	-31.372	-27.301	1.00	16.58	C

ATOM	14838	CB	ILE B 431	-54.400	-31.738	-26.375	1.00	16.53	C
ATOM	14840	CG1	ILE B 431	-54.225	-33.123	-25.728	1.00	16.06	C
ATOM	14843	CD1	ILE B 431	-54.069	-34.212	-26.698	1.00	16.05	C
ATOM	14847	CG2	ILE B 431	-55.701	-31.666	-27.149	1.00	15.49	C
ATOM	14851	C	ILE B 431	-53.431	-29.949	-27.748	1.00	16.35	C
ATOM	14852	O	ILE B 431	-53.856	-29.712	-28.860	1.00	16.74	O
ATOM	14854	N	PHE B 432	-53.110	-29.003	-26.881	1.00	16.15	N
ATOM	14855	CA	PHE B 432	-53.353	-27.577	-27.145	1.00	16.19	C
ATOM	14857	CB	PHE B 432	-52.811	-26.776	-25.956	1.00	16.36	C
ATOM	14860	CG	PHE B 432	-53.007	-25.295	-26.043	1.00	16.54	C
ATOM	14861	CD1	PHE B 432	-53.869	-24.708	-26.946	1.00	16.55	C
ATOM	14863	CE1	PHE B 432	-54.007	-23.334	-26.974	1.00	17.89	C
ATOM	14865	CZ	PHE B 432	-53.305	-22.538	-26.080	1.00	17.96	C
ATOM	14867	CE2	PHE B 432	-52.461	-23.114	-25.169	1.00	17.30	C
ATOM	14869	CD2	PHE B 432	-52.323	-24.482	-25.151	1.00	17.54	C
ATOM	14871	C	PHE B 432	-52.726	-27.104	-28.452	1.00	15.93	C
ATOM	14872	O	PHE B 432	-53.398	-26.525	-29.303	1.00	15.63	O
ATOM	14874	N	ARG B 433	-51.438	-27.387	-28.594	1.00	15.79	N
ATOM	14875	CA	ARG B 433	-50.674	-27.044	-29.778	1.00	15.74	C
ATOM	14877	CB	ARG B 433	-49.196	-27.279	-29.470	1.00	15.68	C
ATOM	14880	CG	ARG B 433	-48.259	-27.218	-30.665	1.00	16.39	C
ATOM	14883	CD	ARG B 433	-48.348	-25.919	-31.388	1.00	16.27	C
ATOM	14886	NE	ARG B 433	-47.969	-24.850	-30.492	1.00	17.34	N
ATOM	14888	CZ	ARG B 433	-48.210	-23.565	-30.703	1.00	17.89	C
ATOM	14889	NH1	ARG B 433	-48.834	-23.142	-31.803	1.00	16.79	N
ATOM	14892	NH2	ARG B 433	-47.805	-22.701	-29.792	1.00	18.53	N
ATOM	14895	C	ARG B 433	-51.114	-27.852	-31.011	1.00	15.85	C
ATOM	14896	O	ARG B 433	-51.197	-27.328	-32.120	1.00	15.40	O
ATOM	14898	N	LEU B 434	-51.381	-29.134	-30.813	1.00	16.18	N
ATOM	14899	CA	LEU B 434	-51.794	-29.986	-31.902	1.00	16.58	C
ATOM	14901	CB	LEU B 434	-51.802	-31.453	-31.468	1.00	16.79	C
ATOM	14904	CG	LEU B 434	-50.408	-32.055	-31.232	1.00	17.49	C
ATOM	14906	CD1	LEU B 434	-50.486	-33.548	-30.864	1.00	16.63	C
ATOM	14910	CD2	LEU B 434	-49.511	-31.836	-32.467	1.00	18.00	C
ATOM	14914	C	LEU B 434	-53.158	-29.559	-32.423	1.00	17.07	C
ATOM	14915	O	LEU B 434	-53.333	-29.442	-33.633	1.00	17.44	O
ATOM	14917	N	CYS B 435	-54.118	-29.308	-31.532	1.00	17.23	N
ATOM	14918	CA	CYS B 435	-55.419	-28.778	-31.947	1.00	17.44	C
ATOM	14920	CB	CYS B 435	-56.301	-28.529	-30.746	1.00	17.38	C
ATOM	14923	SG	CYS B 435	-56.825	-29.993	-29.971	1.00	17.90	S
ATOM	14925	C	CYS B 435	-55.295	-27.460	-32.696	1.00	17.78	C
ATOM	14926	O	CYS B 435	-55.978	-27.232	-33.701	1.00	17.44	O
ATOM	14928	N	ASN B 436	-54.442	-26.585	-32.178	1.00	18.31	N
ATOM	14929	CA	ASN B 436	-54.248	-25.267	-32.769	1.00	19.14	C
ATOM	14931	CB	ASN B 436	-53.392	-24.394	-31.844	1.00	19.31	C
ATOM	14934	CG	ASN B 436	-52.984	-23.068	-32.478	1.00	20.07	C
ATOM	14935	OD1	ASN B 436	-53.456	-22.686	-33.555	1.00	22.72	O
ATOM	14936	ND2	ASN B 436	-52.099	-22.356	-31.800	1.00	20.90	N
ATOM	14939	C	ASN B 436	-53.615	-25.345	-34.157	1.00	19.51	C
ATOM	14940	O	ASN B 436	-54.173	-24.840	-35.131	1.00	19.79	O
ATOM	14942	N	ASP B 437	-52.446	-25.965	-34.245	1.00	19.85	N
ATOM	14943	CA	ASP B 437	-51.754	-26.058	-35.517	1.00	20.01	C
ATOM	14945	CB	ASP B 437	-50.340	-26.638	-35.331	1.00	20.00	C
ATOM	14948	CG	ASP B 437	-49.397	-25.675	-34.557	1.00	21.23	C
ATOM	14949	OD1	ASP B 437	-49.902	-24.749	-33.877	1.00	22.90	O
ATOM	14950	OD2	ASP B 437	-48.149	-25.821	-34.632	1.00	22.18	O
ATOM	14951	C	ASP B 437	-52.629	-26.838	-36.512	1.00	20.00	C
ATOM	14952	O	ASP B 437	-52.678	-26.491	-37.695	1.00	20.02	O
ATOM	14954	N	LEU B 438	-53.372	-27.836	-36.021	1.00	20.02	N
ATOM	14955	CA	LEU B 438	-54.352	-28.558	-36.858	1.00	20.18	C

ATOM	14957	CB	LEU B 438	-55.144	-29.593	-36.049	1.00	19.88	C
ATOM	14960	CG	LEU B 438	-54.699	-31.050	-36.142	1.00	18.97	C
ATOM	14962	CD1	LEU B 438	-55.537	-31.869	-35.200	1.00	18.86	C
ATOM	14966	CD2	LEU B 438	-54.818	-31.585	-37.550	1.00	16.23	C
ATOM	14970	C	LEU B 438	-55.358	-27.644	-37.564	1.00	20.60	C
ATOM	14971	O	LEU B 438	-55.717	-27.889	-38.712	1.00	20.30	O
ATOM	14973	N	ALA B 439	-55.825	-26.620	-36.853	1.00	21.39	N
ATOM	14974	CA	ALA B 439	-56.802	-25.654	-37.377	1.00	21.96	C
ATOM	14976	CB	ALA B 439	-57.306	-24.741	-36.242	1.00	21.77	C
ATOM	14980	C	ALA B 439	-56.214	-24.805	-38.494	1.00	22.47	C
ATOM	14981	O	ALA B 439	-56.850	-24.559	-39.517	1.00	22.13	O
ATOM	14983	N	SER B 440	-54.983	-24.370	-38.282	1.00	23.40	N
ATOM	14984	CA	SER B 440	-54.337	-23.428	-39.169	1.00	24.40	C
ATOM	14986	CB	SER B 440	-53.423	-22.529	-38.342	1.00	24.44	C
ATOM	14989	OG	SER B 440	-52.994	-23.206	-37.166	1.00	25.19	O
ATOM	14991	C	SER B 440	-53.544	-24.103	-40.289	1.00	25.15	C
ATOM	14992	O	SER B 440	-53.093	-23.424	-41.210	1.00	25.24	O
ATOM	14994	N	ALA B 441	-53.392	-25.428	-40.222	1.00	25.99	N
ATOM	14995	CA	ALA B 441	-52.488	-26.160	-41.119	1.00	26.53	C
ATOM	14997	CB	ALA B 441	-52.532	-27.655	-40.830	1.00	26.26	C
ATOM	15001	C	ALA B 441	-52.726	-25.903	-42.608	1.00	27.31	C
ATOM	15002	O	ALA B 441	-51.811	-25.489	-43.310	1.00	27.26	O
ATOM	15004	N	SER B 442	-53.940	-26.134	-43.098	1.00	28.45	N
ATOM	15005	CA	SER B 442	-54.138	-26.138	-44.545	1.00	29.47	C
ATOM	15007	CB	SER B 442	-55.477	-26.769	-44.936	1.00	29.40	C
ATOM	15010	OG	SER B 442	-56.454	-25.784	-45.181	1.00	30.36	O
ATOM	15012	C	SER B 442	-53.966	-24.731	-45.125	1.00	30.22	C
ATOM	15013	O	SER B 442	-53.304	-24.562	-46.139	1.00	30.26	O
ATOM	15015	N	ALA B 443	-54.537	-23.729	-44.467	1.00	31.38	N
ATOM	15016	CA	ALA B 443	-54.308	-22.333	-44.835	1.00	32.22	C
ATOM	15018	CB	ALA B 443	-55.057	-21.410	-43.893	1.00	32.28	C
ATOM	15022	C	ALA B 443	-52.820	-21.992	-44.818	1.00	33.12	C
ATOM	15023	O	ALA B 443	-52.302	-21.419	-45.766	1.00	33.24	O
ATOM	15025	N	GLU B 444	-52.136	-22.351	-43.736	1.00	34.43	N
ATOM	15026	CA	GLU B 444	-50.711	-22.033	-43.582	1.00	35.35	C
ATOM	15028	CB	GLU B 444	-50.237	-22.271	-42.134	1.00	35.46	C
ATOM	15031	CG	GLU B 444	-50.757	-21.227	-41.126	1.00	36.42	C
ATOM	15034	CD	GLU B 444	-50.178	-21.376	-39.708	1.00	37.82	C
ATOM	15035	OE1	GLU B 444	-49.965	-22.506	-39.224	1.00	37.92	O
ATOM	15036	OE2	GLU B 444	-49.955	-20.344	-39.049	1.00	40.09	O
ATOM	15037	C	GLU B 444	-49.840	-22.794	-44.594	1.00	35.87	C
ATOM	15038	O	GLU B 444	-48.870	-22.239	-45.107	1.00	35.85	O
ATOM	15040	N	ILE B 445	-50.193	-24.048	-44.883	1.00	36.63	N
ATOM	15041	CA	ILE B 445	-49.517	-24.837	-45.928	1.00	37.21	C
ATOM	15043	CB	ILE B 445	-49.854	-26.366	-45.832	1.00	37.15	C
ATOM	15045	CG1	ILE B 445	-49.181	-26.991	-44.609	1.00	36.85	C
ATOM	15048	CD1	ILE B 445	-49.867	-28.244	-44.102	1.00	36.62	C
ATOM	15052	CG2	ILE B 445	-49.416	-27.120	-47.093	1.00	36.66	C
ATOM	15056	C	ILE B 445	-49.933	-24.289	-47.290	1.00	37.88	C
ATOM	15057	O	ILE B 445	-50.926	-24.733	-47.874	1.00	38.20	O
ATOM	15059	N	ALA B 446	-49.175	-23.314	-47.785	1.00	38.56	N
ATOM	15060	CA	ALA B 446	-49.544	-22.575	-48.992	1.00	39.03	C
ATOM	15062	CB	ALA B 446	-49.833	-23.522	-50.169	1.00	39.19	C
ATOM	15066	C	ALA B 446	-50.768	-21.740	-48.685	1.00	39.42	C
ATOM	15067	O	ALA B 446	-51.866	-22.288	-48.742	1.00	39.26	O
ATOM	15069	N	ARG B 447	-50.650	-20.443	-48.361	1.00	40.06	N
ATOM	15070	CA	ARG B 447	-49.426	-19.589	-48.361	1.00	40.54	C
ATOM	15072	CB	ARG B 447	-49.276	-18.901	-46.986	1.00	40.70	C
ATOM	15075	CG	ARG B 447	-50.121	-17.634	-46.836	1.00	41.82	C
ATOM	15078	CD	ARG B 447	-50.292	-17.234	-45.377	1.00	42.99	C

ATOM	15081	NE	ARG B 447	-51.628	-17.546	-44.864	1.00	44.22	N
ATOM	15083	CZ	ARG B 447	-51.924	-17.766	-43.580	1.00	46.00	C
ATOM	15084	NH1	ARG B 447	-50.977	-17.730	-42.633	1.00	46.34	N
ATOM	15087	NH2	ARG B 447	-53.184	-18.039	-43.236	1.00	46.39	N
ATOM	15090	C	ARG B 447	-48.095	-20.208	-48.837	1.00	40.35	C
ATOM	15091	O	ARG B 447	-47.897	-20.403	-50.038	1.00	40.95	O
ATOM	15093	N	GLY B 448	-47.170	-20.442	-47.918	1.00	39.87	N
ATOM	15094	CA	GLY B 448	-46.020	-21.306	-48.162	1.00	39.52	C
ATOM	15097	C	GLY B 448	-45.258	-21.517	-46.865	1.00	39.26	C
ATOM	15098	O	GLY B 448	-44.071	-21.842	-46.885	1.00	39.15	O
ATOM	15100	N	GLU B 449	-45.966	-21.344	-45.742	1.00	38.85	N
ATOM	15101	CA	GLU B 449	-45.365	-21.215	-44.426	1.00	38.45	C
ATOM	15103	CB	GLU B 449	-46.288	-20.472	-43.444	1.00	38.76	C
ATOM	15106	CG	GLU B 449	-46.346	-18.943	-43.636	1.00	40.34	C
ATOM	15109	CD	GLU B 449	-47.486	-18.257	-42.839	1.00	42.09	C
ATOM	15110	OE1	GLU B 449	-47.897	-18.787	-41.780	1.00	43.23	O
ATOM	15111	OE2	GLU B 449	-47.972	-17.184	-43.274	1.00	42.13	O
ATOM	15112	C	GLU B 449	-45.073	-22.600	-43.906	1.00	37.54	C
ATOM	15113	O	GLU B 449	-45.810	-23.553	-44.167	1.00	37.34	O
ATOM	15115	N	THR B 450	-43.986	-22.688	-43.159	1.00	36.47	N
ATOM	15116	CA	THR B 450	-43.459	-23.950	-42.702	1.00	35.65	C
ATOM	15118	CB	THR B 450	-42.017	-24.105	-43.245	1.00	35.81	C
ATOM	15120	OG1	THR B 450	-41.593	-25.465	-43.109	1.00	37.67	O
ATOM	15122	CG2	THR B 450	-41.029	-23.152	-42.538	1.00	36.21	C
ATOM	15126	C	THR B 450	-43.528	-24.084	-41.166	1.00	34.22	C
ATOM	15127	O	THR B 450	-43.165	-25.127	-40.619	1.00	34.00	O
ATOM	15129	N	ALA B 451	-44.024	-23.034	-40.495	1.00	32.69	N
ATOM	15130	CA	ALA B 451	-44.144	-22.974	-39.026	1.00	31.28	C
ATOM	15132	CB	ALA B 451	-43.952	-21.546	-38.545	1.00	31.06	C
ATOM	15136	C	ALA B 451	-45.495	-23.500	-38.547	1.00	29.97	C
ATOM	15137	O	ALA B 451	-46.337	-22.735	-38.082	1.00	30.08	O
ATOM	15139	N	ASN B 452	-45.690	-24.809	-38.656	1.00	28.19	N
ATOM	15140	CA	ASN B 452	-46.935	-25.444	-38.272	1.00	26.73	C
ATOM	15142	CB	ASN B 452	-47.929	-25.339	-39.422	1.00	26.45	C
ATOM	15145	CG	ASN B 452	-49.311	-25.832	-39.057	1.00	25.38	C
ATOM	15146	OD1	ASN B 452	-49.608	-27.004	-39.195	1.00	24.57	O
ATOM	15147	ND2	ASN B 452	-50.171	-24.928	-38.618	1.00	24.44	N
ATOM	15150	C	ASN B 452	-46.631	-26.889	-37.952	1.00	26.05	C
ATOM	15151	O	ASN B 452	-45.930	-27.544	-38.695	1.00	26.21	O
ATOM	15153	N	SER B 453	-47.131	-27.388	-36.834	1.00	25.30	N
ATOM	15154	CA	SER B 453	-46.843	-28.760	-36.432	1.00	24.72	C
ATOM	15156	CB	SER B 453	-47.638	-29.139	-35.174	1.00	24.81	C
ATOM	15159	OG	SER B 453	-47.143	-28.455	-34.031	1.00	24.36	O
ATOM	15161	C	SER B 453	-47.111	-29.770	-37.543	1.00	24.12	C
ATOM	15162	O	SER B 453	-46.325	-30.679	-37.743	1.00	24.33	O
ATOM	15164	N	VAL B 454	-48.202	-29.606	-38.272	1.00	23.57	N
ATOM	15165	CA	VAL B 454	-48.546	-30.541	-39.328	1.00	23.27	C
ATOM	15167	CB	VAL B 454	-50.001	-30.361	-39.795	1.00	23.06	C
ATOM	15169	CG1	VAL B 454	-50.363	-31.385	-40.829	1.00	22.41	C
ATOM	15173	CG2	VAL B 454	-50.946	-30.489	-38.622	1.00	22.80	C
ATOM	15177	C	VAL B 454	-47.581	-30.473	-40.521	1.00	23.83	C
ATOM	15178	O	VAL B 454	-47.370	-31.493	-41.175	1.00	24.48	O
ATOM	15180	N	SER B 455	-46.981	-29.309	-40.803	1.00	23.99	N
ATOM	15181	CA	SER B 455	-45.978	-29.198	-41.881	1.00	24.22	C
ATOM	15183	CB	SER B 455	-45.577	-27.767	-42.111	1.00	24.02	C
ATOM	15186	OG	SER B 455	-46.655	-27.098	-42.689	1.00	25.57	O
ATOM	15188	C	SER B 455	-44.708	-29.949	-41.582	1.00	24.65	C
ATOM	15189	O	SER B 455	-44.211	-30.700	-42.417	1.00	25.23	O
ATOM	15191	N	CYS B 456	-44.151	-29.700	-40.405	1.00	24.82	N
ATOM	15192	CA	CYS B 456	-43.024	-30.462	-39.935	1.00	24.81	C

ATOM	15194	CB	CYS B 456	-42.765	-30.170	-38.470	1.00	24.92	C
ATOM	15197	SG	CYS B 456	-42.139	-28.541	-38.190	1.00	25.18	S
ATOM	15199	C	CYS B 456	-43.304	-31.938	-40.124	1.00	24.84	C
ATOM	15200	O	CYS B 456	-42.507	-32.634	-40.728	1.00	25.34	O
ATOM	15202	N	TYR B 457	-44.442	-32.421	-39.644	1.00	24.75	N
ATOM	15203	CA	TYR B 457	-44.735	-33.838	-39.778	1.00	24.85	C
ATOM	15205	CB	TYR B 457	-46.069	-34.197	-39.142	1.00	24.64	C
ATOM	15208	CG	TYR B 457	-46.057	-35.537	-38.451	1.00	24.00	C
ATOM	15209	CD1	TYR B 457	-45.846	-35.635	-37.083	1.00	24.19	C
ATOM	15211	CE1	TYR B 457	-45.848	-36.873	-36.436	1.00	23.72	C
ATOM	15213	CZ	TYR B 457	-46.059	-38.025	-37.170	1.00	23.27	C
ATOM	15214	OH	TYR B 457	-46.063	-39.244	-36.554	1.00	22.94	O
ATOM	15216	CE2	TYR B 457	-46.268	-37.961	-38.528	1.00	23.68	C
ATOM	15218	CD2	TYR B 457	-46.265	-36.711	-39.163	1.00	24.16	C
ATOM	15220	C	TYR B 457	-44.676	-34.267	-41.251	1.00	25.30	C
ATOM	15221	O	TYR B 457	-43.965	-35.208	-41.584	1.00	24.94	O
ATOM	15223	N	MET B 458	-45.388	-33.569	-42.134	1.00	26.02	N
ATOM	15224	CA	MET B 458	-45.207	-33.775	-43.585	1.00	26.76	C
ATOM	15226	CB	MET B 458	-45.913	-32.685	-44.407	1.00	26.84	C
ATOM	15229	CG	MET B 458	-47.419	-32.720	-44.429	1.00	27.32	C
ATOM	15232	SD	MET B 458	-48.092	-31.238	-45.229	1.00	27.76	S
ATOM	15233	CE	MET B 458	-47.283	-31.311	-46.830	1.00	26.95	C
ATOM	15237	C	MET B 458	-43.723	-33.759	-43.996	1.00	27.06	C
ATOM	15238	O	MET B 458	-43.192	-34.746	-44.485	1.00	26.94	O
ATOM	15240	N	ARG B 459	-43.069	-32.624	-43.800	1.00	27.58	N
ATOM	15241	CA	ARG B 459	-41.720	-32.419	-44.289	1.00	28.38	C
ATOM	15243	CB	ARG B 459	-41.313	-30.965	-44.041	1.00	28.96	C
ATOM	15246	CG	ARG B 459	-39.862	-30.642	-44.328	1.00	31.93	C
ATOM	15249	CD	ARG B 459	-39.656	-29.139	-44.480	1.00	36.88	C
ATOM	15252	NE	ARG B 459	-40.318	-28.349	-43.427	1.00	41.70	N
ATOM	15254	CZ	ARG B 459	-39.890	-28.241	-42.161	1.00	45.68	C
ATOM	15255	NH1	ARG B 459	-38.795	-28.899	-41.751	1.00	48.08	N
ATOM	15258	NH2	ARG B 459	-40.566	-27.489	-41.284	1.00	45.92	N
ATOM	15261	C	ARG B 459	-40.717	-33.401	-43.672	1.00	28.03	C
ATOM	15262	O	ARG B 459	-39.829	-33.881	-44.356	1.00	28.13	O
ATOM	15264	N	THR B 460	-40.882	-33.711	-42.393	1.00	27.88	N
ATOM	15265	CA	THR B 460	-40.002	-34.635	-41.664	1.00	27.75	C
ATOM	15267	CB	THR B 460	-40.240	-34.502	-40.134	1.00	27.79	C
ATOM	15269	OG1	THR B 460	-39.805	-33.209	-39.689	1.00	27.98	O
ATOM	15271	CG2	THR B 460	-39.519	-35.589	-39.356	1.00	27.46	C
ATOM	15275	C	THR B 460	-40.171	-36.114	-42.045	1.00	27.72	C
ATOM	15276	O	THR B 460	-39.211	-36.857	-42.069	1.00	27.60	O
ATOM	15278	N	LYS B 461	-41.398	-36.538	-42.313	1.00	28.05	N
ATOM	15279	CA	LYS B 461	-41.701	-37.929	-42.649	1.00	28.25	C
ATOM	15281	CB	LYS B 461	-42.972	-38.381	-41.914	1.00	28.32	C
ATOM	15284	CG	LYS B 461	-42.767	-38.763	-40.453	1.00	28.57	C
ATOM	15287	CD	LYS B 461	-42.569	-40.271	-40.299	1.00	29.79	C
ATOM	15290	CE	LYS B 461	-41.913	-40.651	-38.970	1.00	30.32	C
ATOM	15293	NZ	LYS B 461	-42.651	-40.184	-37.749	1.00	30.60	N
ATOM	15297	C	LYS B 461	-41.876	-38.139	-44.155	1.00	28.50	C
ATOM	15298	O	LYS B 461	-42.071	-39.269	-44.598	1.00	28.58	O
ATOM	15300	N	GLY B 462	-41.818	-37.057	-44.933	1.00	28.79	N
ATOM	15301	CA	GLY B 462	-41.978	-37.114	-46.385	1.00	29.04	C
ATOM	15304	C	GLY B 462	-43.333	-37.610	-46.853	1.00	29.39	C
ATOM	15305	O	GLY B 462	-43.410	-38.452	-47.738	1.00	29.53	O
ATOM	15307	N	ILE B 463	-44.406	-37.081	-46.270	1.00	29.91	N
ATOM	15308	CA	ILE B 463	-45.764	-37.552	-46.572	1.00	30.27	C
ATOM	15310	CB	ILE B 463	-46.333	-38.444	-45.433	1.00	30.13	C
ATOM	15312	CG1	ILE B 463	-46.308	-37.706	-44.092	1.00	30.22	C
ATOM	15315	CD1	ILE B 463	-46.983	-38.462	-42.965	1.00	30.14	C

ATOM	15319	CG2 ILE B 463	-45.558	-39.740	-45.340	1.00	29.85	C
ATOM	15323	C ILE B 463	-46.775	-36.429	-46.881	1.00	30.75	C
ATOM	15324	O ILE B 463	-46.595	-35.270	-46.495	1.00	30.52	O
ATOM	15326	N SER B 464	-47.842	-36.822	-47.581	1.00	31.33	N
ATOM	15327	CA SER B 464	-48.951	-35.944	-47.952	1.00	31.68	C
ATOM	15329	CB SER B 464	-49.948	-36.709	-48.830	1.00	32.02	C
ATOM	15332	OG SER B 464	-50.547	-37.793	-48.123	1.00	32.84	O
ATOM	15334	C SER B 464	-49.705	-35.396	-46.747	1.00	31.62	C
ATOM	15335	O SER B 464	-49.897	-36.097	-45.747	1.00	31.69	O
ATOM	15337	N GLU B 465	-50.165	-34.153	-46.874	1.00	31.42	N
ATOM	15338	CA GLU B 465	-50.881	-33.466	-45.803	1.00	31.20	C
ATOM	15340	CB GLU B 465	-51.465	-32.141	-46.324	1.00	31.33	C
ATOM	15343	CG GLU B 465	-52.441	-31.455	-45.359	1.00	31.42	C
ATOM	15346	CD GLU B 465	-52.768	-30.006	-45.712	1.00	31.30	C
ATOM	15347	OE1 GLU B 465	-52.325	-29.490	-46.765	1.00	30.56	O
ATOM	15348	OE2 GLU B 465	-53.487	-29.380	-44.906	1.00	31.84	O
ATOM	15349	C GLU B 465	-51.985	-34.325	-45.202	1.00	30.97	C
ATOM	15350	O GLU B 465	-52.175	-34.327	-43.992	1.00	30.72	O
ATOM	15352	N GLU B 466	-52.709	-35.052	-46.055	1.00	30.98	N
ATOM	15353	CA GLU B 466	-53.865	-35.833	-45.618	1.00	30.65	C
ATOM	15355	CB GLU B 466	-54.553	-36.541	-46.788	1.00	30.62	C
ATOM	15358	CG GLU B 466	-55.733	-37.401	-46.336	1.00	31.50	C
ATOM	15361	CD GLU B 466	-56.353	-38.240	-47.446	1.00	33.14	C
ATOM	15362	OE1 GLU B 466	-56.512	-37.736	-48.581	1.00	34.40	O
ATOM	15363	OE2 GLU B 466	-56.699	-39.413	-47.173	1.00	33.21	O
ATOM	15364	C GLU B 466	-53.477	-36.851	-44.564	1.00	30.17	C
ATOM	15365	O GLU B 466	-54.215	-37.049	-43.606	1.00	30.46	O
ATOM	15367	N LEU B 467	-52.333	-37.506	-44.726	1.00	29.59	N
ATOM	15368	CA LEU B 467	-51.953	-38.540	-43.769	1.00	29.19	C
ATOM	15370	CB LEU B 467	-51.640	-39.869	-44.476	1.00	29.27	C
ATOM	15373	CG LEU B 467	-50.356	-40.081	-45.268	1.00	29.84	C
ATOM	15375	CD1 LEU B 467	-49.341	-40.761	-44.357	1.00	31.06	C
ATOM	15379	CD2 LEU B 467	-50.589	-40.920	-46.523	1.00	29.47	C
ATOM	15383	C LEU B 467	-50.865	-38.092	-42.794	1.00	28.54	C
ATOM	15384	O LEU B 467	-50.490	-38.834	-41.894	1.00	28.28	O
ATOM	15386	N ALA B 468	-50.388	-36.862	-42.962	1.00	28.04	N
ATOM	15387	CA ALA B 468	-49.682	-36.158	-41.893	1.00	27.54	C
ATOM	15389	CB ALA B 468	-48.922	-34.970	-42.435	1.00	27.34	C
ATOM	15393	C ALA B 468	-50.711	-35.698	-40.865	1.00	27.17	C
ATOM	15394	O ALA B 468	-50.482	-35.801	-39.662	1.00	27.40	O
ATOM	15396	N THR B 469	-51.841	-35.182	-41.349	1.00	26.57	N
ATOM	15397	CA THR B 469	-52.957	-34.795	-40.490	1.00	25.87	C
ATOM	15399	CB THR B 469	-54.130	-34.297	-41.328	1.00	25.53	C
ATOM	15401	OG1 THR B 469	-53.737	-33.111	-42.025	1.00	24.99	O
ATOM	15403	CG2 THR B 469	-55.316	-33.993	-40.453	1.00	24.62	C
ATOM	15407	C THR B 469	-53.412	-35.988	-39.669	1.00	25.95	C
ATOM	15408	O THR B 469	-53.551	-35.925	-38.450	1.00	25.80	O
ATOM	15410	N GLU B 470	-53.618	-37.091	-40.367	1.00	26.06	N
ATOM	15411	CA GLU B 470	-54.043	-38.344	-39.756	1.00	26.10	C
ATOM	15413	CB GLU B 470	-54.057	-39.436	-40.848	1.00	26.60	C
ATOM	15416	CG GLU B 470	-55.038	-40.587	-40.646	1.00	28.01	C
ATOM	15419	CD GLU B 470	-54.553	-41.897	-41.286	1.00	29.05	C
ATOM	15420	OE1 GLU B 470	-54.298	-41.930	-42.515	1.00	27.96	O
ATOM	15421	OE2 GLU B 470	-54.435	-42.892	-40.536	1.00	30.42	O
ATOM	15422	C GLU B 470	-53.148	-38.746	-38.557	1.00	25.28	C
ATOM	15423	O GLU B 470	-53.629	-39.266	-37.556	1.00	25.08	O
ATOM	15425	N SER B 471	-51.850	-38.497	-38.654	1.00	24.70	N
ATOM	15426	CA SER B 471	-50.930	-38.906	-37.599	1.00	24.26	C
ATOM	15428	CB SER B 471	-49.494	-38.842	-38.085	1.00	24.11	C
ATOM	15431	OG SER B 471	-49.315	-39.693	-39.186	1.00	23.98	O

ATOM	15433	C	SER B 471	-51.085	-38.029	-36.377	1.00	23.92	C
ATOM	15434	O	SER B 471	-51.240	-38.530	-35.270	1.00	24.09	O
ATOM	15436	N	VAL B 472	-51.032	-36.719	-36.582	1.00	23.48	N
ATOM	15437	CA	VAL B 472	-51.279	-35.774	-35.510	1.00	23.15	C
ATOM	15439	CB	VAL B 472	-51.377	-34.332	-36.045	1.00	23.03	C
ATOM	15441	CG1	VAL B 472	-51.739	-33.347	-34.929	1.00	21.74	C
ATOM	15445	CG2	VAL B 472	-50.067	-33.945	-36.728	1.00	22.27	C
ATOM	15449	C	VAL B 472	-52.572	-36.192	-34.830	1.00	23.49	C
ATOM	15450	O	VAL B 472	-52.663	-36.276	-33.614	1.00	23.29	O
ATOM	15452	N	MET B 473	-53.573	-36.508	-35.626	1.00	24.34	N
ATOM	15453	CA	MET B 473	-54.850	-36.901	-35.070	1.00	24.93	C
ATOM	15455	CB	MET B 473	-55.816	-37.266	-36.190	1.00	25.09	C
ATOM	15458	CG	MET B 473	-57.191	-36.788	-35.912	1.00	26.94	C
ATOM	15461	SD	MET B 473	-57.328	-35.035	-36.251	1.00	29.69	S
ATOM	15462	CE	MET B 473	-58.316	-35.080	-37.775	1.00	29.12	C
ATOM	15466	C	MET B 473	-54.668	-38.072	-34.100	1.00	25.00	C
ATOM	15467	O	MET B 473	-55.135	-38.018	-32.963	1.00	24.78	O
ATOM	15469	N	ASN B 474	-53.965	-39.110	-34.556	1.00	25.30	N
ATOM	15470	CA	ASN B 474	-53.681	-40.293	-33.738	1.00	25.55	C
ATOM	15472	CB	ASN B 474	-53.143	-41.439	-34.601	1.00	25.82	C
ATOM	15475	CG	ASN B 474	-54.243	-42.137	-35.402	1.00	27.62	C
ATOM	15476	OD1	ASN B 474	-54.238	-42.102	-36.638	1.00	29.33	O
ATOM	15477	ND2	ASN B 474	-55.196	-42.780	-34.697	1.00	28.58	N
ATOM	15480	C	ASN B 474	-52.711	-40.034	-32.592	1.00	25.20	C
ATOM	15481	O	ASN B 474	-52.711	-40.773	-31.621	1.00	25.71	O
ATOM	15483	N	LEU B 475	-51.883	-38.999	-32.704	1.00	24.65	N
ATOM	15484	CA	LEU B 475	-50.987	-38.595	-31.614	1.00	24.04	C
ATOM	15486	CB	LEU B 475	-49.939	-37.626	-32.142	1.00	23.88	C
ATOM	15489	CG	LEU B 475	-48.837	-37.285	-31.164	1.00	23.70	C
ATOM	15491	CD1	LEU B 475	-47.936	-38.471	-31.004	1.00	23.97	C
ATOM	15495	CD2	LEU B 475	-48.088	-36.104	-31.679	1.00	23.99	C
ATOM	15499	C	LEU B 475	-51.747	-37.940	-30.453	1.00	23.65	C
ATOM	15500	O	LEU B 475	-51.323	-38.014	-29.299	1.00	23.74	O
ATOM	15502	N	ILE B 476	-52.858	-37.282	-30.766	1.00	23.15	N
ATOM	15503	CA	ILE B 476	-53.728	-36.711	-29.742	1.00	22.50	C
ATOM	15505	CB	ILE B 476	-54.779	-35.755	-30.359	1.00	22.27	C
ATOM	15507	CG1	ILE B 476	-54.099	-34.508	-30.910	1.00	20.55	C
ATOM	15510	CD1	ILE B 476	-55.000	-33.671	-31.693	1.00	19.19	C
ATOM	15514	CG2	ILE B 476	-55.817	-35.363	-29.336	1.00	22.15	C
ATOM	15518	C	ILE B 476	-54.398	-37.849	-28.985	1.00	22.40	C
ATOM	15519	O	ILE B 476	-54.316	-37.900	-27.767	1.00	22.32	O
ATOM	15521	N	ASP B 477	-55.023	-38.777	-29.708	1.00	22.37	N
ATOM	15522	CA	ASP B 477	-55.616	-39.975	-29.090	1.00	22.55	C
ATOM	15524	CB	ASP B 477	-55.996	-41.007	-30.151	1.00	22.67	C
ATOM	15527	CG	ASP B 477	-57.262	-40.636	-30.892	1.00	24.93	C
ATOM	15528	OD1	ASP B 477	-57.690	-39.458	-30.818	1.00	28.87	O
ATOM	15529	OD2	ASP B 477	-57.847	-41.523	-31.549	1.00	27.92	O
ATOM	15530	C	ASP B 477	-54.675	-40.622	-28.083	1.00	22.10	C
ATOM	15531	O	ASP B 477	-55.048	-40.877	-26.938	1.00	21.83	O
ATOM	15533	N	GLU B 478	-53.444	-40.867	-28.510	1.00	21.73	N
ATOM	15534	CA	GLU B 478	-52.472	-41.506	-27.645	1.00	21.69	C
ATOM	15536	CB	GLU B 478	-51.233	-41.920	-28.437	1.00	21.97	C
ATOM	15539	CG	GLU B 478	-50.786	-43.360	-28.150	1.00	24.18	C
ATOM	15542	CD	GLU B 478	-51.793	-44.429	-28.618	1.00	27.03	C
ATOM	15543	OE1	GLU B 478	-51.599	-45.627	-28.285	1.00	29.15	O
ATOM	15544	OE2	GLU B 478	-52.770	-44.080	-29.323	1.00	28.37	O
ATOM	15545	C	GLU B 478	-52.107	-40.608	-26.449	1.00	20.97	C
ATOM	15546	O	GLU B 478	-51.928	-41.110	-25.316	1.00	21.15	O
ATOM	15548	N	THR B 479	-52.031	-39.291	-26.687	1.00	19.82	N
ATOM	15549	CA	THR B 479	-51.787	-38.319	-25.604	1.00	18.59	C

ATOM	15551	CB	THR B 479	-51.587	-36.897	-26.129	1.00	17.95	C
ATOM	15553	OG1	THR B 479	-50.291	-36.779	-26.712	1.00	17.43	O
ATOM	15555	CG2	THR B 479	-51.654	-35.938	-25.009	1.00	17.90	C
ATOM	15559	C	THR B 479	-52.927	-38.342	-24.580	1.00	18.16	C
ATOM	15560	O	THR B 479	-52.695	-38.339	-23.383	1.00	18.06	O
ATOM	15562	N	TRP B 480	-54.160	-38.390	-25.057	1.00	17.77	N
ATOM	15563	CA	TRP B 480	-55.301	-38.493	-24.169	1.00	17.54	C
ATOM	15565	CB	TRP B 480	-56.605	-38.464	-24.970	1.00	17.49	C
ATOM	15568	CG	TRP B 480	-57.239	-37.099	-25.028	1.00	17.53	C
ATOM	15569	CD1	TRP B 480	-57.150	-36.186	-26.041	1.00	17.04	C
ATOM	15571	NE1	TRP B 480	-57.857	-35.067	-25.721	1.00	16.57	N
ATOM	15573	CE2	TRP B 480	-58.423	-35.236	-24.484	1.00	16.41	C
ATOM	15574	CD2	TRP B 480	-58.055	-36.498	-24.020	1.00	16.68	C
ATOM	15575	CE3	TRP B 480	-58.507	-36.914	-22.765	1.00	16.92	C
ATOM	15577	CZ3	TRP B 480	-59.300	-36.076	-22.039	1.00	16.97	C
ATOM	15579	CH2	TRP B 480	-59.654	-34.823	-22.524	1.00	16.91	C
ATOM	15581	CZ2	TRP B 480	-59.223	-34.383	-23.742	1.00	16.89	C
ATOM	15583	C	TRP B 480	-55.232	-39.760	-23.334	1.00	17.75	C
ATOM	15584	O	TRP B 480	-55.582	-39.747	-22.164	1.00	17.74	O
ATOM	15586	N	LYS B 481	-54.784	-40.861	-23.933	1.00	18.09	N
ATOM	15587	CA	LYS B 481	-54.721	-42.124	-23.210	1.00	18.10	C
ATOM	15589	CB	LYS B 481	-54.277	-43.282	-24.115	1.00	18.01	C
ATOM	15592	CG	LYS B 481	-55.311	-43.803	-25.114	1.00	16.89	C
ATOM	15595	CD	LYS B 481	-54.613	-44.663	-26.199	1.00	15.83	C
ATOM	15598	CE	LYS B 481	-55.587	-45.500	-27.025	1.00	14.74	C
ATOM	15601	NZ	LYS B 481	-55.119	-45.725	-28.413	1.00	12.34	N
ATOM	15605	C	LYS B 481	-53.750	-41.970	-22.054	1.00	18.46	C
ATOM	15606	O	LYS B 481	-53.976	-42.494	-20.969	1.00	18.44	O
ATOM	15608	N	LYS B 482	-52.662	-41.252	-22.279	1.00	18.80	N
ATOM	15609	CA	LYS B 482	-51.727	-41.019	-21.192	1.00	19.58	C
ATOM	15611	CB	LYS B 482	-50.425	-40.424	-21.727	1.00	20.03	C
ATOM	15614	CG	LYS B 482	-49.499	-41.500	-22.277	1.00	22.03	C
ATOM	15617	CD	LYS B 482	-48.706	-41.063	-23.506	1.00	24.87	C
ATOM	15620	CE	LYS B 482	-47.908	-42.268	-24.066	1.00	26.65	C
ATOM	15623	NZ	LYS B 482	-47.287	-41.996	-25.408	1.00	28.80	N
ATOM	15627	C	LYS B 482	-52.363	-40.156	-20.096	1.00	19.56	C
ATOM	15628	O	LYS B 482	-52.367	-40.537	-18.926	1.00	19.19	O
ATOM	15630	N	MET B 483	-52.922	-39.013	-20.483	1.00	20.00	N
ATOM	15631	CA	MET B 483	-53.677	-38.165	-19.552	1.00	20.37	C
ATOM	15633	CB	MET B 483	-54.426	-37.057	-20.295	1.00	20.25	C
ATOM	15636	CG	MET B 483	-53.529	-35.951	-20.807	1.00	20.52	C
ATOM	15639	SD	MET B 483	-54.461	-34.483	-21.225	1.00	20.47	S
ATOM	15640	CE	MET B 483	-55.361	-35.084	-22.649	1.00	22.07	C
ATOM	15644	C	MET B 483	-54.668	-38.981	-18.741	1.00	20.72	C
ATOM	15645	O	MET B 483	-54.809	-38.779	-17.543	1.00	20.37	O
ATOM	15647	N	ASN B 484	-55.345	-39.910	-19.403	1.00	21.57	N
ATOM	15648	CA	ASN B 484	-56.339	-40.731	-18.744	1.00	22.46	C
ATOM	15650	CB	ASN B 484	-57.042	-41.631	-19.765	1.00	22.23	C
ATOM	15653	CG	ASN B 484	-58.118	-40.900	-20.576	1.00	21.61	C
ATOM	15654	OD1	ASN B 484	-58.412	-39.722	-20.362	1.00	20.25	O
ATOM	15655	ND2	ASN B 484	-58.720	-41.622	-21.510	1.00	20.80	N
ATOM	15658	C	ASN B 484	-55.723	-41.555	-17.598	1.00	23.83	C
ATOM	15659	O	ASN B 484	-56.262	-41.591	-16.484	1.00	23.76	O
ATOM	15661	N	LYS B 485	-54.585	-42.190	-17.868	1.00	25.69	N
ATOM	15662	CA	LYS B 485	-53.854	-42.937	-16.843	1.00	27.11	C
ATOM	15664	CB	LYS B 485	-52.679	-43.712	-17.447	1.00	27.13	C
ATOM	15667	CG	LYS B 485	-52.073	-44.755	-16.489	1.00	28.85	C
ATOM	15670	CD	LYS B 485	-51.643	-46.082	-17.187	1.00	31.17	C
ATOM	15673	CE	LYS B 485	-50.150	-46.132	-17.555	1.00	32.11	C
ATOM	15676	NZ	LYS B 485	-49.291	-46.282	-16.344	1.00	33.08	N

ATOM	15680	C	LYS B 485	-53.370	-42.035	-15.701	1.00	28.37	C
ATOM	15681	O	LYS B 485	-53.288	-42.481	-14.563	1.00	28.59	O
ATOM	15683	N	GLU B 486	-53.083	-40.766	-15.980	1.00	29.89	N
ATOM	15684	CA	GLU B 486	-52.683	-39.849	-14.910	1.00	31.09	C
ATOM	15686	CB	GLU B 486	-52.174	-38.511	-15.470	1.00	31.55	C
ATOM	15689	CG	GLU B 486	-51.059	-37.851	-14.643	1.00	32.69	C
ATOM	15692	CD	GLU B 486	-49.785	-38.689	-14.605	1.00	34.51	C
ATOM	15693	OE1	GLU B 486	-49.104	-38.682	-13.557	1.00	36.35	O
ATOM	15694	OE2	GLU B 486	-49.476	-39.368	-15.611	1.00	34.55	O
ATOM	15695	C	GLU B 486	-53.836	-39.615	-13.939	1.00	31.70	C
ATOM	15696	O	GLU B 486	-53.684	-39.868	-12.750	1.00	31.90	O
ATOM	15698	N	LYS B 487	-54.983	-39.164	-14.453	1.00	32.71	N
ATOM	15699	CA	LYS B 487	-56.177	-38.856	-13.629	1.00	33.57	C
ATOM	15701	CB	LYS B 487	-57.343	-38.402	-14.519	1.00	33.54	C
ATOM	15704	CG	LYS B 487	-58.703	-38.196	-13.823	1.00	33.45	C
ATOM	15707	CD	LYS B 487	-58.686	-37.024	-12.866	1.00	33.71	C
ATOM	15710	CE	LYS B 487	-59.906	-36.996	-11.931	1.00	34.48	C
ATOM	15713	NZ	LYS B 487	-61.123	-36.378	-12.526	1.00	35.03	N
ATOM	15717	C	LYS B 487	-56.626	-40.046	-12.802	1.00	34.60	C
ATOM	15718	O	LYS B 487	-57.126	-39.892	-11.681	1.00	34.55	O
ATOM	15720	N	LEU B 488	-56.453	-41.233	-13.371	1.00	35.94	N
ATOM	15721	CA	LEU B 488	-56.878	-42.451	-12.723	1.00	37.00	C
ATOM	15723	CB	LEU B 488	-57.280	-43.489	-13.778	1.00	36.90	C
ATOM	15726	CG	LEU B 488	-58.081	-44.714	-13.312	1.00	37.52	C
ATOM	15728	CD1	LEU B 488	-58.758	-44.536	-11.946	1.00	38.76	C
ATOM	15732	CD2	LEU B 488	-59.121	-45.084	-14.346	1.00	38.42	C
ATOM	15736	C	LEU B 488	-55.786	-42.980	-11.803	1.00	38.02	C
ATOM	15737	O	LEU B 488	-56.063	-43.315	-10.653	1.00	38.19	O
ATOM	15739	N	GLY B 489	-54.550	-43.016	-12.298	1.00	39.44	N
ATOM	15740	CA	GLY B 489	-53.441	-43.710	-11.621	1.00	40.62	C
ATOM	15743	C	GLY B 489	-52.563	-42.837	-10.740	1.00	41.76	C
ATOM	15744	O	GLY B 489	-51.357	-42.703	-10.988	1.00	41.83	O
ATOM	15746	N	GLY B 490	-53.183	-42.230	-9.726	1.00	43.10	N
ATOM	15747	CA	GLY B 490	-52.481	-41.538	-8.636	1.00	43.82	C
ATOM	15750	C	GLY B 490	-51.131	-40.924	-8.967	1.00	44.31	C
ATOM	15751	O	GLY B 490	-50.084	-41.560	-8.831	1.00	44.25	O
ATOM	15753	N	SER B 491	-51.161	-39.672	-9.396	1.00	44.87	N
ATOM	15754	CA	SER B 491	-49.939	-38.892	-9.575	1.00	45.22	C
ATOM	15756	CB	SER B 491	-50.224	-37.723	-10.537	1.00	45.34	C
ATOM	15759	OG	SER B 491	-51.411	-37.032	-10.166	1.00	45.38	O
ATOM	15761	C	SER B 491	-49.419	-38.389	-8.203	1.00	45.14	C
ATOM	15762	O	SER B 491	-49.896	-38.824	-7.140	1.00	45.11	O
ATOM	15764	N	LEU B 492	-48.429	-37.491	-8.227	1.00	44.77	N
ATOM	15765	CA	LEU B 492	-48.046	-36.752	-7.018	1.00	44.31	C
ATOM	15767	CB	LEU B 492	-46.726	-35.971	-7.188	1.00	44.55	C
ATOM	15770	CG	LEU B 492	-45.530	-36.546	-7.972	1.00	45.73	C
ATOM	15772	CD1	LEU B 492	-44.373	-35.556	-7.859	1.00	46.25	C
ATOM	15776	CD2	LEU B 492	-45.084	-37.972	-7.530	1.00	46.14	C
ATOM	15780	C	LEU B 492	-49.151	-35.760	-6.703	1.00	43.30	C
ATOM	15781	O	LEU B 492	-49.315	-35.370	-5.555	1.00	43.49	O
ATOM	15783	N	PHE B 493	-49.903	-35.367	-7.735	1.00	42.05	N
ATOM	15784	CA	PHE B 493	-50.822	-34.239	-7.668	1.00	41.03	C
ATOM	15786	CB	PHE B 493	-50.913	-33.549	-9.028	1.00	40.81	C
ATOM	15789	CG	PHE B 493	-49.696	-32.771	-9.408	1.00	38.96	C
ATOM	15790	CD1	PHE B 493	-49.651	-31.408	-9.233	1.00	37.06	C
ATOM	15792	CE1	PHE B 493	-48.543	-30.697	-9.603	1.00	36.41	C
ATOM	15794	CZ	PHE B 493	-47.471	-31.341	-10.158	1.00	35.87	C
ATOM	15796	CE2	PHE B 493	-47.507	-32.690	-10.343	1.00	36.39	C
ATOM	15798	CD2	PHE B 493	-48.612	-33.400	-9.977	1.00	37.36	C
ATOM	15800	C	PHE B 493	-52.240	-34.603	-7.268	1.00	40.76	C

ATOM 15801	O	PHE B 493	-52.686	-35.733	-7.419	1.00	40.35	O
ATOM 15803	N	ALA B 494	-52.943	-33.584	-6.790	1.00	40.72	N
ATOM 15804	CA	ALA B 494	-54.355	-33.658	-6.471	1.00	40.68	C
ATOM 15806	CB	ALA B 494	-54.786	-32.362	-5.770	1.00	40.70	C
ATOM 15810	C	ALA B 494	-55.175	-33.851	-7.744	1.00	40.59	C
ATOM 15811	O	ALA B 494	-55.013	-33.099	-8.716	1.00	40.95	O
ATOM 15813	N	LYS B 495	-56.071	-34.835	-7.733	1.00	40.12	N
ATOM 15814	CA	LYS B 495	-57.000	-35.047	-8.854	1.00	39.67	C
ATOM 15816	CB	LYS B 495	-58.008	-36.166	-8.520	1.00	39.98	C
ATOM 15819	CG	LYS B 495	-57.388	-37.590	-8.496	1.00	40.57	C
ATOM 15822	CD	LYS B 495	-58.468	-38.688	-8.419	1.00	41.18	C
ATOM 15825	CE	LYS B 495	-57.878	-40.065	-8.132	1.00	41.28	C
ATOM 15828	NZ	LYS B 495	-58.736	-41.159	-8.671	1.00	41.39	N
ATOM 15832	C	LYS B 495	-57.717	-33.752	-9.333	1.00	38.82	C
ATOM 15833	O	LYS B 495	-57.847	-33.535	-10.537	1.00	38.70	O
ATOM 15835	N	PRO B 496	-58.169	-32.890	-8.397	1.00	37.71	N
ATOM 15836	CA	PRO B 496	-58.701	-31.551	-8.687	1.00	36.80	C
ATOM 15838	CB	PRO B 496	-58.820	-30.937	-7.304	1.00	36.94	C
ATOM 15841	CG	PRO B 496	-59.255	-32.073	-6.483	1.00	37.90	C
ATOM 15844	CD	PRO B 496	-58.534	-33.294	-7.030	1.00	37.86	C
ATOM 15847	C	PRO B 496	-57.845	-30.625	-9.536	1.00	35.57	C
ATOM 15848	O	PRO B 496	-58.389	-29.693	-10.139	1.00	35.93	O
ATOM 15849	N	PHE B 497	-56.527	-30.836	-9.546	1.00	33.52	N
ATOM 15850	CA	PHE B 497	-55.654	-30.083	-10.439	1.00	31.34	C
ATOM 15852	CB	PHE B 497	-54.355	-29.710	-9.750	1.00	30.86	C
ATOM 15855	CG	PHE B 497	-53.394	-28.989	-10.634	1.00	29.04	C
ATOM 15856	CD1	PHE B 497	-53.639	-27.700	-11.022	1.00	27.89	C
ATOM 15858	CE1	PHE B 497	-52.752	-27.028	-11.834	1.00	27.39	C
ATOM 15860	CZ	PHE B 497	-51.599	-27.653	-12.259	1.00	27.40	C
ATOM 15862	CE2	PHE B 497	-51.341	-28.941	-11.877	1.00	27.44	C
ATOM 15864	CD2	PHE B 497	-52.236	-29.604	-11.074	1.00	28.18	C
ATOM 15866	C	PHE B 497	-55.390	-30.870	-11.714	1.00	30.24	C
ATOM 15867	O	PHE B 497	-55.301	-30.266	-12.786	1.00	30.39	O
ATOM 15869	N	VAL B 498	-55.295	-32.200	-11.630	1.00	28.50	N
ATOM 15870	CA	VAL B 498	-55.076	-32.983	-12.854	1.00	27.54	C
ATOM 15872	CB	VAL B 498	-54.778	-34.472	-12.611	1.00	27.34	C
ATOM 15874	CG1	VAL B 498	-56.027	-35.194	-12.262	1.00	28.14	C
ATOM 15878	CG2	VAL B 498	-53.731	-34.649	-11.530	1.00	27.04	C
ATOM 15882	C	VAL B 498	-56.277	-32.857	-13.784	1.00	26.60	C
ATOM 15883	O	VAL B 498	-56.129	-32.864	-14.995	1.00	26.33	O
ATOM 15885	N	GLU B 499	-57.468	-32.731	-13.214	1.00	25.87	N
ATOM 15886	CA	GLU B 499	-58.662	-32.484	-14.022	1.00	25.20	C
ATOM 15888	CB	GLU B 499	-59.975	-32.639	-13.206	1.00	25.26	C
ATOM 15891	CG	GLU B 499	-61.287	-32.694	-14.028	1.00	25.16	C
ATOM 15894	CD	GLU B 499	-61.422	-33.922	-14.973	1.00	26.79	C
ATOM 15895	OE1	GLU B 499	-60.607	-34.880	-14.928	1.00	26.39	O
ATOM 15896	OE2	GLU B 499	-62.379	-33.926	-15.785	1.00	27.83	O
ATOM 15897	C	GLU B 499	-58.542	-31.097	-14.641	1.00	24.27	C
ATOM 15898	O	GLU B 499	-58.819	-30.948	-15.832	1.00	24.32	O
ATOM 15900	N	THR B 500	-58.098	-30.093	-13.882	1.00	22.90	N
ATOM 15901	CA	THR B 500	-58.027	-28.764	-14.487	1.00	22.29	C
ATOM 15903	CB	THR B 500	-57.908	-27.592	-13.488	1.00	22.43	C
ATOM 15905	OG1	THR B 500	-56.536	-27.259	-13.274	1.00	22.65	O
ATOM 15907	CG2	THR B 500	-58.614	-27.899	-12.183	1.00	22.84	C
ATOM 15911	C	THR B 500	-56.940	-28.671	-15.565	1.00	21.20	C
ATOM 15912	O	THR B 500	-57.120	-27.939	-16.536	1.00	21.06	O
ATOM 15914	N	ALA B 501	-55.845	-29.414	-15.422	1.00	19.88	N
ATOM 15915	CA	ALA B 501	-54.906	-29.574	-16.537	1.00	19.24	C
ATOM 15917	CB	ALA B 501	-53.704	-30.393	-16.115	1.00	19.22	C
ATOM 15921	C	ALA B 501	-55.581	-30.226	-17.757	1.00	18.68	C

ATOM 15922 O ALA B 501	-55.427 -29.768 -18.881	1.00 18.32	O
ATOM 15924 N ILE B 502	-56.324 -31.300 -17.523	1.00 18.11	N
ATOM 15925 CA ILE B 502	-57.023 -31.988 -18.591	1.00 17.68	C
ATOM 15927 CB ILE B 502	-57.750 -33.251 -18.059	1.00 17.63	C
ATOM 15929 CG1 ILE B 502	-56.713 -34.350 -17.769	1.00 18.16	C
ATOM 15932 CD1 ILE B 502	-57.253 -35.645 -17.115	1.00 17.02	C
ATOM 15936 CG2 ILE B 502	-58.756 -33.762 -19.063	1.00 16.93	C
ATOM 15940 C ILE B 502	-57.990 -31.030 -19.289	1.00 17.50	C
ATOM 15941 O ILE B 502	-58.153 -31.071 -20.512	1.00 17.62	O
ATOM 15943 N ASN B 503	-58.604 -30.139 -18.524	1.00 17.23	N
ATOM 15944 CA ASN B 503	-59.545 -29.184 -19.102	1.00 17.07	C
ATOM 15946 CB ASN B 503	-60.184 -28.327 -18.016	1.00 17.21	C
ATOM 15949 CG ASN B 503	-61.155 -29.094 -17.176	1.00 17.29	C
ATOM 15950 OD1 ASN B 503	-61.721 -30.097 -17.609	1.00 16.66	O
ATOM 15951 ND2 ASN B 503	-61.365 -28.621 -15.958	1.00 19.16	N
ATOM 15954 C ASN B 503	-58.920 -28.270 -20.149	1.00 16.59	C
ATOM 15955 O ASN B 503	-59.611 -27.811 -21.071	1.00 16.67	O
ATOM 15957 N LEU B 504	-57.629 -27.992 -20.014	1.00 15.79	N
ATOM 15958 CA LEU B 504	-56.948 -27.204 -21.031	1.00 15.44	C
ATOM 15960 CB LEU B 504	-55.530 -26.865 -20.609	1.00 15.07	C
ATOM 15963 CG LEU B 504	-54.933 -25.837 -21.552	1.00 14.21	C
ATOM 15965 CD1 LEU B 504	-54.246 -24.726 -20.800	1.00 14.21	C
ATOM 15969 CD2 LEU B 504	-54.000 -26.538 -22.462	1.00 14.27	C
ATOM 15973 C LEU B 504	-56.954 -27.933 -22.375	1.00 15.59	C
ATOM 15974 O LEU B 504	-57.177 -27.323 -23.414	1.00 15.30	O
ATOM 15976 N ALA B 505	-56.732 -29.242 -22.341	1.00 15.89	N
ATOM 15977 CA ALA B 505	-56.914 -30.073 -23.514	1.00 16.23	C
ATOM 15979 CB ALA B 505	-56.608 -31.483 -23.188	1.00 16.05	C
ATOM 15983 C ALA B 505	-58.352 -29.955 -24.010	1.00 16.69	C
ATOM 15984 O ALA B 505	-58.598 -29.738 -25.204	1.00 16.88	O
ATOM 15986 N ARG B 506	-59.300 -30.071 -23.090	1.00 16.88	N
ATOM 15987 CA ARG B 506	-60.709 -29.972 -23.457	1.00 17.20	C
ATOM 15989 CB ARG B 506	-61.630 -30.170 -22.248	1.00 17.34	C
ATOM 15992 CG ARG B 506	-61.549 -31.537 -21.614	1.00 17.46	C
ATOM 15995 CD ARG B 506	-62.837 -31.917 -20.962	1.00 17.48	C
ATOM 15998 NE ARG B 506	-62.783 -33.280 -20.435	1.00 18.37	N
ATOM 16000 CZ ARG B 506	-62.404 -33.613 -19.201	1.00 19.03	C
ATOM 16001 NH1 ARG B 506	-62.008 -32.690 -18.324	1.00 18.83	N
ATOM 16004 NH2 ARG B 506	-62.415 -34.892 -18.841	1.00 19.72	N
ATOM 16007 C ARG B 506	-61.038 -28.646 -24.114	1.00 17.32	C
ATOM 16008 O ARG B 506	-61.754 -28.622 -25.111	1.00 17.41	O
ATOM 16010 N GLN B 507	-60.540 -27.543 -23.564	1.00 17.59	N
ATOM 16011 CA GLN B 507	-60.864 -26.233 -24.139	1.00 18.05	C
ATOM 16013 CB GLN B 507	-60.382 -25.090 -23.258	1.00 18.16	C
ATOM 16016 CG GLN B 507	-60.798 -23.709 -23.764	1.00 17.57	C
ATOM 16019 CD GLN B 507	-62.291 -23.508 -23.723	1.00 17.42	C
ATOM 16020 OE1 GLN B 507	-62.953 -23.959 -22.800	1.00 17.55	O
ATOM 16021 NE2 GLN B 507	-62.832 -22.823 -24.722	1.00 17.56	N
ATOM 16024 C GLN B 507	-60.258 -26.060 -25.521	1.00 18.49	C
ATOM 16025 O GLN B 507	-60.855 -25.413 -26.380	1.00 18.73	O
ATOM 16027 N SER B 508	-59.066 -26.619 -25.724	1.00 18.85	N
ATOM 16028 CA SER B 508	-58.427 -26.621 -27.037	1.00 19.07	C
ATOM 16030 CB SER B 508	-57.108 -27.364 -26.969	1.00 19.04	C
ATOM 16033 OG SER B 508	-56.304 -26.772 -25.979	1.00 20.10	O
ATOM 16035 C SER B 508	-59.305 -27.296 -28.065	1.00 19.12	C
ATOM 16036 O SER B 508	-59.438 -26.827 -29.187	1.00 19.00	O
ATOM 16038 N HIS B 509	-59.905 -28.409 -27.674	1.00 19.26	N
ATOM 16039 CA HIS B 509	-60.814 -29.095 -28.559	1.00 19.33	C
ATOM 16041 CB HIS B 509	-61.275 -30.412 -27.959	1.00 19.30	C
ATOM 16044 CG HIS B 509	-60.263 -31.501 -28.067	1.00 18.71	C

ATOM	16045	ND1 HIS B 509	-59.941 -32.089 -29.267	1.00 18.57	N
ATOM	16047	CE1 HIS B 509	-59.020 -33.012 -29.064	1.00 19.07	C
ATOM	16049	NE2 HIS B 509	-58.738 -33.044 -27.774	1.00 18.52	N
ATOM	16051	CD2 HIS B 509	-59.499 -32.105 -27.129	1.00 18.43	C
ATOM	16053	C HIS B 509	-62.017 -28.251 -28.878	1.00 19.60	C
ATOM	16054	O HIS B 509	-62.491 -28.286 -29.969	1.00 19.68	O
ATOM	16056	N CYS B 510	-62.540 -27.502 -27.932	1.00 20.12	N
ATOM	16057	CA CYS B 510	-63.748 -26.738 -28.216	1.00 20.63	C
ATOM	16059	CB CYS B 510	-64.493 -26.492 -26.915	1.00 20.70	C
ATOM	16062	SG CYS B 510	-64.856 -28.020 -26.065	1.00 22.93	S
ATOM	16064	C CYS B 510	-63.465 -25.422 -28.950	1.00 20.55	C
ATOM	16065	O CYS B 510	-64.338 -24.881 -29.622	1.00 20.16	O
ATOM	16067	N THR B 511	-62.236 -24.933 -28.827	1.00 20.89	N
ATOM	16068	CA THR B 511	-61.833 -23.654 -29.378	1.00 21.32	C
ATOM	16070	CB THR B 511	-60.682 -23.107 -28.550	1.00 20.76	C
ATOM	16072	OG1 THR B 511	-61.208 -22.619 -27.324	1.00 19.50	O
ATOM	16074	CG2 THR B 511	-59.961 -21.994 -29.253	1.00 19.79	C
ATOM	16078	C THR B 511	-61.417 -23.735 -30.851	1.00 23.11	C
ATOM	16079	O THR B 511	-61.910 -22.992 -31.694	1.00 22.75	O
ATOM	16081	N TYR B 512	-60.501 -24.646 -31.152	1.00 25.46	N
ATOM	16082	CA TYR B 512	-59.867 -24.704 -32.462	1.00 27.09	C
ATOM	16084	CB TYR B 512	-58.397 -25.107 -32.327	1.00 27.16	C
ATOM	16087	CG TYR B 512	-57.598 -24.056 -31.582	1.00 27.26	C
ATOM	16088	CD1 TYR B 512	-57.300 -22.834 -32.182	1.00 28.18	C
ATOM	16090	CE1 TYR B 512	-56.583 -21.844 -31.511	1.00 27.98	C
ATOM	16092	CZ TYR B 512	-56.164 -22.068 -30.226	1.00 27.72	C
ATOM	16093	OH TYR B 512	-55.462 -21.072 -29.584	1.00 27.09	O
ATOM	16095	CE2 TYR B 512	-56.453 -23.281 -29.602	1.00 27.52	C
ATOM	16097	CD2 TYR B 512	-57.171 -24.262 -30.278	1.00 26.68	C
ATOM	16099	C TYR B 512	-60.633 -25.623 -33.382	1.00 28.79	C
ATOM	16100	O TYR B 512	-61.050 -25.179 -34.444	1.00 28.98	O
ATOM	16102	N HIS B 513	-60.789 -26.893 -32.991	1.00 31.06	N
ATOM	16103	CA HIS B 513	-61.865 -27.800 -33.492	1.00 33.07	C
ATOM	16105	CB HIS B 513	-62.740 -28.208 -32.265	1.00 33.66	C
ATOM	16108	CG HIS B 513	-64.181 -28.602 -32.535	1.00 34.73	C
ATOM	16109	ND1 HIS B 513	-65.035 -27.903 -33.367	1.00 35.70	N
ATOM	16111	CE1 HIS B 513	-66.234 -28.464 -33.343	1.00 35.40	C
ATOM	16113	NE2 HIS B 513	-66.206 -29.469 -32.491	1.00 35.10	N
ATOM	16115	CD2 HIS B 513	-64.943 -29.565 -31.956	1.00 35.07	C
ATOM	16117	C HIS B 513	-62.692 -27.179 -34.623	1.00 34.03	C
ATOM	16118	O HIS B 513	-63.003 -27.848 -35.627	1.00 34.55	O
ATOM	16120	N ASN B 514	-63.047 -25.906 -34.444	1.00 34.65	N
ATOM	16121	CA ASN B 514	-63.713 -25.123 -35.483	1.00 35.15	C
ATOM	16123	CB ASN B 514	-63.845 -23.643 -35.070	1.00 35.01	C
ATOM	16126	CG ASN B 514	-64.704 -23.456 -33.820	1.00 33.48	C
ATOM	16127	OD1 ASN B 514	-64.787 -24.343 -32.968	1.00 30.92	O
ATOM	16128	ND2 ASN B 514	-65.327 -22.295 -33.703	1.00 32.16	N
ATOM	16131	C ASN B 514	-63.192 -25.274 -36.942	1.00 36.08	C
ATOM	16132	O ASN B 514	-62.010 -25.051 -37.254	1.00 35.69	O
ATOM	16134	N GLY B 515	-64.119 -25.765 -37.772	1.00 37.15	N
ATOM	16135	CA GLY B 515	-64.224 -25.493 -39.177	1.00 37.94	C
ATOM	16138	C GLY B 515	-65.269 -24.384 -39.284	1.00 39.01	C
ATOM	16139	O GLY B 515	-64.904 -23.286 -39.702	1.00 39.58	O
ATOM	16141	N ASP B 516	-66.550 -24.585 -38.909	1.00 39.95	N
ATOM	16142	CA ASP B 516	-67.167 -25.810 -38.349	1.00 40.65	C
ATOM	16144	CB ASP B 516	-67.296 -25.652 -36.836	1.00 41.09	C
ATOM	16147	CG ASP B 516	-66.318 -26.511 -36.063	1.00 43.63	C
ATOM	16148	OD1 ASP B 516	-65.697 -27.434 -36.662	1.00 45.83	O
ATOM	16149	OD2 ASP B 516	-66.152 -26.246 -34.843	1.00 47.12	O
ATOM	16150	C ASP B 516	-68.603 -26.067 -38.853	1.00 40.78	C

ATOM	16151	O	ASP B 516	-69.134	-25.301	-39.660	1.00	41.28	O
ATOM	16153	N	ALA B 517	-69.222	-27.141	-38.349	1.00	40.66	N
ATOM	16154	CA	ALA B 517	-70.675	-27.396	-38.468	1.00	40.60	C
ATOM	16156	CB	ALA B 517	-71.422	-26.622	-37.359	1.00	40.27	C
ATOM	16160	C	ALA B 517	-71.295	-27.106	-39.862	1.00	40.79	C
ATOM	16161	O	ALA B 517	-70.654	-27.296	-40.899	1.00	40.67	O
ATOM	16163	N	HIS B 518	-72.562	-26.696	-39.874	1.00	41.13	N
ATOM	16164	CA	HIS B 518	-73.162	-26.033	-41.039	1.00	41.60	C
ATOM	16166	CB	HIS B 518	-74.446	-26.756	-41.480	1.00	42.24	C
ATOM	16169	CG	HIS B 518	-74.207	-28.179	-41.910	1.00	45.54	C
ATOM	16170	ND1	HIS B 518	-73.441	-28.509	-43.012	1.00	48.51	N
ATOM	16172	CE1	HIS B 518	-73.392	-29.825	-43.140	1.00	48.79	C
ATOM	16174	NE2	HIS B 518	-74.092	-30.364	-42.156	1.00	49.12	N
ATOM	16176	CD2	HIS B 518	-74.608	-29.357	-41.369	1.00	48.28	C
ATOM	16178	C	HIS B 518	-73.377	-24.545	-40.686	1.00	40.84	C
ATOM	16179	O	HIS B 518	-74.487	-23.997	-40.769	1.00	40.60	O
ATOM	16181	N	THR B 519	-72.261	-23.932	-40.275	1.00	40.10	N
ATOM	16182	CA	THR B 519	-72.143	-22.521	-39.878	1.00	39.39	C
ATOM	16184	CB	THR B 519	-72.561	-22.263	-38.393	1.00	39.40	C
ATOM	16186	OG1	THR B 519	-71.823	-23.131	-37.517	1.00	39.55	O
ATOM	16188	CG2	THR B 519	-74.074	-22.460	-38.178	1.00	38.98	C
ATOM	16192	C	THR B 519	-70.654	-22.164	-40.059	1.00	38.81	C
ATOM	16193	O	THR B 519	-69.800	-23.050	-40.067	1.00	38.48	O
ATOM	16195	N	SER B 520	-70.338	-20.881	-40.199	1.00	38.19	N
ATOM	16196	CA	SER B 520	-68.959	-20.453	-40.521	1.00	37.77	C
ATOM	16198	CB	SER B 520	-68.983	-18.976	-40.960	1.00	37.76	C
ATOM	16201	OG	SER B 520	-68.760	-18.106	-39.870	1.00	38.49	O
ATOM	16203	C	SER B 520	-67.974	-20.724	-39.340	1.00	37.24	C
ATOM	16204	O	SER B 520	-68.394	-21.267	-38.320	1.00	36.89	O
ATOM	16206	N	PRO B 521	-66.671	-20.350	-39.472	1.00	36.96	N
ATOM	16207	CA	PRO B 521	-65.692	-20.638	-38.396	1.00	36.80	C
ATOM	16209	CB	PRO B 521	-64.320	-20.354	-39.041	1.00	36.68	C
ATOM	16212	CG	PRO B 521	-64.599	-19.593	-40.309	1.00	37.19	C
ATOM	16215	CD	PRO B 521	-66.091	-19.469	-40.507	1.00	37.10	C
ATOM	16218	C	PRO B 521	-65.915	-19.763	-37.168	1.00	36.68	C
ATOM	16219	O	PRO B 521	-66.233	-20.289	-36.105	1.00	36.90	O
ATOM	16220	N	ASP B 522	-65.744	-18.445	-37.298	1.00	36.61	N
ATOM	16221	CA	ASP B 522	-66.373	-17.528	-36.347	1.00	36.58	C
ATOM	16223	CB	ASP B 522	-65.938	-16.070	-36.543	1.00	36.82	C
ATOM	16226	CG	ASP B 522	-64.456	-15.846	-36.243	1.00	37.94	C
ATOM	16227	OD1	ASP B 522	-63.848	-16.627	-35.473	1.00	39.62	O
ATOM	16228	OD2	ASP B 522	-63.891	-14.876	-36.790	1.00	39.40	O
ATOM	16229	C	ASP B 522	-67.836	-17.718	-36.681	1.00	36.12	C
ATOM	16230	O	ASP B 522	-68.150	-18.306	-37.702	1.00	35.98	O
ATOM	16232	N	GLU B 523	-68.732	-17.249	-35.830	1.00	35.69	N
ATOM	16233	CA	GLU B 523	-70.152	-17.605	-35.927	1.00	35.46	C
ATOM	16235	CB	GLU B 523	-70.703	-17.560	-37.379	1.00	35.51	C
ATOM	16238	CG	GLU B 523	-70.535	-16.175	-38.061	1.00	36.11	C
ATOM	16241	CD	GLU B 523	-70.854	-16.136	-39.572	1.00	36.72	C
ATOM	16242	OE1	GLU B 523	-71.848	-16.767	-40.013	1.00	37.20	O
ATOM	16243	OE2	GLU B 523	-70.101	-15.452	-40.314	1.00	35.53	O
ATOM	16244	C	GLU B 523	-70.440	-18.943	-35.216	1.00	34.97	C
ATOM	16245	O	GLU B 523	-71.557	-19.148	-34.765	1.00	35.02	O
ATOM	16247	N	LEU B 524	-69.462	-19.846	-35.099	1.00	34.59	N
ATOM	16248	CA	LEU B 524	-69.503	-20.832	-34.004	1.00	34.45	C
ATOM	16250	CB	LEU B 524	-68.735	-22.136	-34.286	1.00	34.21	C
ATOM	16253	CG	LEU B 524	-69.530	-23.434	-34.508	1.00	33.92	C
ATOM	16255	CD1	LEU B 524	-68.657	-24.627	-34.227	1.00	32.21	C
ATOM	16259	CD2	LEU B 524	-70.794	-23.520	-33.641	1.00	34.06	C
ATOM	16263	C	LEU B 524	-68.901	-20.162	-32.791	1.00	34.38	C

ATOM	16264	O	LEU B 524	-69.535	-20.067	-31.741	1.00	34.49	O
ATOM	16266	N	THR B 525	-67.670	-19.688	-32.958	1.00	34.30	N
ATOM	16267	CA	THR B 525	-66.913	-19.073	-31.874	1.00	34.23	C
ATOM	16269	CB	THR B 525	-65.570	-18.524	-32.380	1.00	34.04	C
ATOM	16271	OG1	THR B 525	-64.894	-19.549	-33.112	1.00	33.83	O
ATOM	16273	CG2	THR B 525	-64.689	-18.087	-31.227	1.00	33.37	C
ATOM	16277	C	THR B 525	-67.702	-17.960	-31.184	1.00	34.46	C
ATOM	16278	O	THR B 525	-67.684	-17.858	-29.955	1.00	34.57	O
ATOM	16280	N	ARG B 526	-68.401	-17.144	-31.969	1.00	34.54	N
ATOM	16281	CA	ARG B 526	-69.200	-16.056	-31.418	1.00	34.60	C
ATOM	16283	CB	ARG B 526	-69.571	-15.059	-32.513	1.00	35.04	C
ATOM	16286	CG	ARG B 526	-70.256	-13.800	-32.013	1.00	37.13	C
ATOM	16289	CD	ARG B 526	-70.143	-12.656	-33.026	1.00	39.87	C
ATOM	16292	NE	ARG B 526	-70.353	-13.075	-34.420	1.00	42.28	N
ATOM	16294	CZ	ARG B 526	-71.543	-13.292	-35.000	1.00	44.26	C
ATOM	16295	NH1	ARG B 526	-72.685	-13.159	-34.320	1.00	44.66	N
ATOM	16298	NH2	ARG B 526	-71.593	-13.658	-36.279	1.00	44.49	N
ATOM	16301	C	ARG B 526	-70.448	-16.593	-30.735	1.00	33.91	C
ATOM	16302	O	ARG B 526	-70.848	-16.074	-29.704	1.00	33.74	O
ATOM	16304	N	LYS B 527	-71.048	-17.637	-31.306	1.00	33.39	N
ATOM	16305	CA	LYS B 527	-72.240	-18.252	-30.722	1.00	33.15	C
ATOM	16307	CB	LYS B 527	-72.837	-19.332	-31.639	1.00	33.42	C
ATOM	16310	CG	LYS B 527	-73.898	-18.804	-32.609	1.00	34.55	C
ATOM	16313	CD	LYS B 527	-74.643	-19.919	-33.358	1.00	35.43	C
ATOM	16316	CE	LYS B 527	-75.230	-19.383	-34.667	1.00	36.11	C
ATOM	16319	NZ	LYS B 527	-76.076	-20.372	-35.384	1.00	37.04	N
ATOM	16323	C	LYS B 527	-71.925	-18.865	-29.378	1.00	32.36	C
ATOM	16324	O	LYS B 527	-72.589	-18.584	-28.387	1.00	32.56	O
ATOM	16326	N	ARG B 528	-70.909	-19.713	-29.358	1.00	31.53	N
ATOM	16327	CA	ARG B 528	-70.502	-20.401	-28.138	1.00	30.77	C
ATOM	16329	CB	ARG B 528	-69.283	-21.286	-28.414	1.00	30.55	C
ATOM	16332	CG	ARG B 528	-69.624	-22.519	-29.252	1.00	29.48	C
ATOM	16335	CD	ARG B 528	-68.418	-23.417	-29.493	1.00	28.05	C
ATOM	16338	NE	ARG B 528	-68.811	-24.774	-29.880	1.00	26.55	N
ATOM	16340	CZ	ARG B 528	-67.968	-25.719	-30.288	1.00	25.87	C
ATOM	16341	NH1	ARG B 528	-66.668	-25.472	-30.385	1.00	26.33	N
ATOM	16344	NH2	ARG B 528	-68.424	-26.920	-30.615	1.00	25.49	N
ATOM	16347	C	ARG B 528	-70.225	-19.415	-27.004	1.00	30.42	C
ATOM	16348	O	ARG B 528	-70.721	-19.582	-25.885	1.00	30.30	O
ATOM	16350	N	VAL B 529	-69.455	-18.376	-27.307	1.00	29.96	N
ATOM	16351	CA	VAL B 529	-69.194	-17.312	-26.342	1.00	29.47	C
ATOM	16353	CB	VAL B 529	-68.261	-16.227	-26.932	1.00	29.40	C
ATOM	16355	CG1	VAL B 529	-68.269	-14.968	-26.081	1.00	29.18	C
ATOM	16359	CG2	VAL B 529	-66.845	-16.782	-27.062	1.00	28.90	C
ATOM	16363	C	VAL B 529	-70.516	-16.723	-25.841	1.00	28.98	C
ATOM	16364	O	VAL B 529	-70.759	-16.682	-24.641	1.00	28.97	O
ATOM	16366	N	LEU B 530	-71.384	-16.308	-26.752	1.00	28.54	N
ATOM	16367	CA	LEU B 530	-72.709	-15.824	-26.353	1.00	28.23	C
ATOM	16369	CB	LEU B 530	-73.631	-15.604	-27.559	1.00	28.09	C
ATOM	16372	CG	LEU B 530	-73.767	-14.150	-28.007	1.00	28.15	C
ATOM	16374	CD1	LEU B 530	-72.412	-13.570	-28.406	1.00	27.99	C
ATOM	16378	CD2	LEU B 530	-74.784	-14.034	-29.150	1.00	28.85	C
ATOM	16382	C	LEU B 530	-73.389	-16.767	-25.367	1.00	27.88	C
ATOM	16383	O	LEU B 530	-73.932	-16.319	-24.363	1.00	28.11	O
ATOM	16385	N	SER B 531	-73.352	-18.066	-25.652	1.00	27.28	N
ATOM	16386	CA	SER B 531	-74.081	-19.049	-24.849	1.00	26.68	C
ATOM	16388	CB	SER B 531	-74.209	-20.363	-25.612	1.00	26.68	C
ATOM	16391	OG	SER B 531	-72.970	-21.033	-25.681	1.00	26.67	O
ATOM	16393	C	SER B 531	-73.410	-19.312	-23.514	1.00	26.20	C
ATOM	16394	O	SER B 531	-74.076	-19.550	-22.511	1.00	25.86	O

ATOM	16396	N	VAL B 532	-72.085	-19.281	-23.512	1.00	25.86	N
ATOM	16397	CA	VAL B 532	-71.320	-19.562	-22.306	1.00	25.51	C
ATOM	16399	CB	VAL B 532	-69.891	-20.055	-22.651	1.00	25.34	C
ATOM	16401	CG1	VAL B 532	-68.955	-19.932	-21.460	1.00	24.01	C
ATOM	16405	CG2	VAL B 532	-69.950	-21.484	-23.158	1.00	24.95	C
ATOM	16409	C	VAL B 532	-71.271	-18.357	-21.373	1.00	25.52	C
ATOM	16410	O	VAL B 532	-71.377	-18.523	-20.164	1.00	25.82	O
ATOM	16412	N	ILE B 533	-71.137	-17.158	-21.935	1.00	25.39	N
ATOM	16413	CA	ILE B 533	-70.875	-15.954	-21.151	1.00	25.48	C
ATOM	16415	CB	ILE B 533	-69.666	-15.185	-21.732	1.00	25.30	C
ATOM	16417	CG1	ILE B 533	-68.375	-15.928	-21.451	1.00	24.91	C
ATOM	16420	CD1	ILE B 533	-68.118	-16.120	-19.989	1.00	25.00	C
ATOM	16424	CG2	ILE B 533	-69.553	-13.797	-21.130	1.00	25.81	C
ATOM	16428	C	ILE B 533	-72.065	-14.983	-21.030	1.00	25.77	C
ATOM	16429	O	ILE B 533	-72.537	-14.707	-19.928	1.00	25.77	O
ATOM	16431	N	THR B 534	-72.541	-14.442	-22.142	1.00	25.98	N
ATOM	16432	CA	THR B 534	-73.446	-13.295	-22.064	1.00	26.31	C
ATOM	16434	CB	THR B 534	-73.094	-12.253	-23.143	1.00	26.21	C
ATOM	16436	OG1	THR B 534	-72.914	-12.909	-24.397	1.00	26.72	O
ATOM	16438	CG2	THR B 534	-71.795	-11.535	-22.781	1.00	25.97	C
ATOM	16442	C	THR B 534	-74.961	-13.621	-22.064	1.00	26.50	C
ATOM	16443	O	THR B 534	-75.713	-12.974	-21.344	1.00	26.67	O
ATOM	16445	N	GLU B 535	-75.410	-14.616	-22.829	1.00	26.55	N
ATOM	16446	CA	GLU B 535	-76.851	-14.892	-22.961	1.00	26.50	C
ATOM	16448	CB	GLU B 535	-77.197	-15.197	-24.426	1.00	26.68	C
ATOM	16451	CG	GLU B 535	-77.226	-13.936	-25.296	1.00	27.68	C
ATOM	16454	CD	GLU B 535	-78.193	-14.014	-26.475	1.00	28.75	C
ATOM	16455	OE1	GLU B 535	-79.370	-13.588	-26.330	1.00	27.72	O
ATOM	16456	OE2	GLU B 535	-77.762	-14.492	-27.549	1.00	30.06	O
ATOM	16457	C	GLU B 535	-77.355	-16.019	-22.056	1.00	26.19	C
ATOM	16458	O	GLU B 535	-77.021	-17.171	-22.279	1.00	26.15	O
ATOM	16460	N	PRO B 536	-78.186	-15.697	-21.047	1.00	26.05	N
ATOM	16461	CA	PRO B 536	-78.692	-16.763	-20.196	1.00	26.05	C
ATOM	16463	CB	PRO B 536	-79.495	-16.018	-19.118	1.00	25.87	C
ATOM	16466	CG	PRO B 536	-79.140	-14.621	-19.239	1.00	25.88	C
ATOM	16469	CD	PRO B 536	-78.760	-14.402	-20.660	1.00	26.21	C
ATOM	16472	C	PRO B 536	-79.608	-17.717	-20.943	1.00	26.14	C
ATOM	16473	O	PRO B 536	-80.173	-17.359	-21.973	1.00	25.99	O
ATOM	16474	N	ILE B 537	-79.746	-18.924	-20.411	1.00	26.38	N
ATOM	16475	CA	ILE B 537	-80.662	-19.905	-20.960	1.00	26.58	C
ATOM	16477	CB	ILE B 537	-80.443	-21.292	-20.333	1.00	26.47	C
ATOM	16479	CG1	ILE B 537	-79.023	-21.789	-20.600	1.00	26.38	C
ATOM	16482	CD1	ILE B 537	-78.709	-23.102	-19.907	1.00	26.40	C
ATOM	16486	CG2	ILE B 537	-81.430	-22.300	-20.890	1.00	26.71	C
ATOM	16490	C	ILE B 537	-82.072	-19.420	-20.657	1.00	26.93	C
ATOM	16491	O	ILE B 537	-82.347	-18.965	-19.545	1.00	26.81	O
ATOM	16493	N	LEU B 538	-82.963	-19.491	-21.641	1.00	27.36	N
ATOM	16494	CA	LEU B 538	-84.309	-18.984	-21.439	1.00	27.68	C
ATOM	16496	CB	LEU B 538	-85.181	-19.094	-22.698	1.00	27.75	C
ATOM	16499	CG	LEU B 538	-84.782	-18.291	-23.949	1.00	27.47	C
ATOM	16501	CD1	LEU B 538	-85.992	-18.109	-24.844	1.00	27.11	C
ATOM	16505	CD2	LEU B 538	-84.164	-16.933	-23.621	1.00	26.99	C
ATOM	16509	C	LEU B 538	-84.919	-19.752	-20.288	1.00	28.01	C
ATOM	16512	N	PRO B 539	-85.645	-19.046	-19.421	1.00	28.69	N
ATOM	16513	CA	PRO B 539	-86.028	-19.631	-18.152	1.00	29.05	C
ATOM	16515	CB	PRO B 539	-86.640	-18.452	-17.397	1.00	29.08	C
ATOM	16518	CG	PRO B 539	-87.144	-17.548	-18.447	1.00	28.86	C
ATOM	16521	CD	PRO B 539	-86.298	-17.748	-19.660	1.00	28.62	C
ATOM	16524	C	PRO B 539	-87.035	-20.767	-18.259	1.00	29.51	C
ATOM	16525	O	PRO B 539	-87.665	-20.970	-19.297	1.00	29.23	O

ATOM	16526	N	PHE B 540	-87.164	-21.500	-17.159	1.00	30.25	N
ATOM	16527	CA	PHE B 540	-88.089	-22.613	-17.070	1.00	30.56	C
ATOM	16529	CB	PHE B 540	-87.971	-23.297	-15.708	1.00	30.73	C
ATOM	16532	CG	PHE B 540	-88.848	-24.499	-15.567	1.00	30.49	C
ATOM	16533	CD1	PHE B 540	-89.910	-24.505	-14.683	1.00	30.32	C
ATOM	16535	CE1	PHE B 540	-90.718	-25.614	-14.568	1.00	30.59	C
ATOM	16537	CZ	PHE B 540	-90.475	-26.722	-15.346	1.00	30.55	C
ATOM	16539	CE2	PHE B 540	-89.420	-26.720	-16.234	1.00	30.68	C
ATOM	16541	CD2	PHE B 540	-88.619	-25.617	-16.342	1.00	30.43	C
ATOM	16543	C	PHE B 540	-89.507	-22.120	-17.257	1.00	30.76	C
ATOM	16544	O	PHE B 540	-89.967	-21.257	-16.508	1.00	30.59	O
ATOM	16546	N	GLU B 541	-90.184	-22.677	-18.259	1.00	31.13	N
ATOM	16547	CA	GLU B 541	-91.553	-22.289	-18.608	1.00	31.45	C
ATOM	16549	CB	GLU B 541	-91.525	-21.186	-19.680	1.00	31.60	C
ATOM	16552	CG	GLU B 541	-92.860	-20.454	-19.911	1.00	32.56	C
ATOM	16555	CD	GLU B 541	-93.773	-21.133	-20.935	1.00	33.64	C
ATOM	16556	OE1	GLU B 541	-93.261	-21.922	-21.771	1.00	35.06	O
ATOM	16557	OE2	GLU B 541	-95.002	-20.868	-20.901	1.00	32.44	O
ATOM	16558	C	GLU B 541	-92.326	-23.510	-19.105	1.00	31.31	C
ATOM	16559	O	GLU B 541	-92.860	-24.286	-18.310	1.00	31.30	O
ATOM	16562	MG	MG C 1	-42.844	11.427	13.309	1.00	46.29	MG
ATOM	16561	MG	MG C 2	-46.615	-18.454	-33.231	1.00	48.45	MG
ATOM	16563	O	HOH E 1	-50.507	-5.408	-4.491	1.00	17.42	O
ATOM	16566	O	HOH E 2	-64.748	-36.007	-11.725	1.00	2.00	O
ATOM	16569	O	HOH E 3	-40.643	-2.220	-34.996	1.00	2.00	O
ATOM	16572	O	HOH E 4	-36.090	-9.757	-37.074	1.00	17.20	O
ATOM	16575	O	HOH E 5	-46.117	-37.662	-22.916	1.00	17.03	O
ATOM	16578	O	HOH E 6	-49.541	35.476	6.921	1.00	14.03	O
ATOM	16581	O	HOH E 7	-32.288	27.572	16.443	1.00	13.16	O
ATOM	16584	O	HOH E 8	-50.706	10.207	-15.061	1.00	18.73	O
ATOM	16587	O	HOH E 9	-77.188	36.767	-9.218	1.00	8.05	O
ATOM	16590	O	HOH E 10	-90.260	-31.248	-21.071	1.00	2.00	O
ATOM	16593	O	HOH E 11	-70.920	-33.414	-3.884	1.00	9.35	O
ATOM	16596	O	HOH E 12	-37.761	-21.294	-9.249	1.00	25.78	O
ATOM	16599	O	HOH E 13	-76.050	23.855	-18.293	1.00	2.00	O
ATOM	16602	O	HOH E 14	-76.876	-19.575	-22.856	1.00	19.25	O
ATOM	16605	O	HOH E 15	-40.936	11.629	-24.832	1.00	22.29	O
ATOM	16608	O	HOH E 16	-85.551	34.832	3.260	1.00	25.14	O
ATOM	16611	O	HOH E 17	-56.825	31.771	-7.464	1.00	17.27	O
ATOM	16614	O	HOH E 18	-76.222	39.261	-6.613	1.00	25.47	O

2009240505 12 Jul 2013

[0478] Various modifications and variations of the described method and system of the invention will be apparent to those skilled in the art without departing from the scope and spirit of the invention. Although the invention has been described in connection with specific preferred embodiments, it should be understood that the invention as claimed should not be unduly limited to such specific embodiments. Indeed, various modifications of the described modes for carrying out the invention, which are obvious to those skilled in the relevant fields, are intended to be within the scope of the following claims.

[0479] Where the terms "comprise", "comprises", "comprised" or "comprising" are used in this specification, they are to be interpreted as specifying the presence of the stated features, integers, steps or components referred to, but not to preclude the presence or addition of one or more other feature, integer, step, component or group thereof.

[0480] Further, any prior art reference or statement provided in the specification is not to be taken as an admission that such art constitutes, or is to be understood as constituting, part of the common general knowledge in Australia.

2009240505 02 Aug 2013

The claims defining the invention are as follows:

1. An isolated isoprene synthase variant comprising one or more amino acid substitution(s) at one or more amino acid positions corresponding to a poplar isoprene synthase having the sequence of SEQ ID NO:120 wherein said substitution(s) are selected from the group consisting of V10M, F12S, T15A, E18G, V58I, V58F, L70Q, L70R, L70V, L70T, T71P, V79L, E89D, G94A, S119F, F120L, G127R, E175V, T212I, S257A, R262G, A266G, F280L, N297K, F305L, L319M, E323K, A328T, D342E, A359T, K366N, E368D, L374M, S396T, V418S, K438N, H440R, T442A, I449V, A469S, K500R, K505Q, G507S, S509N, F511Y, and N532K; and wherein the variant is capable of more effectively converting dimethylallyl diphosphate (DMAPP) to isoprene, as compared to an isoprene synthase variant without a substitution.
2. The isolated isoprene synthase variant of claim 1, wherein the variant further comprises one or more amino acid substitutions selected from the group consisting of G127R/F511Y, L70Q/G94A/R262G/F305L, F12S/T15A/E18G/N297K, S396T/T442I, V10M/E323K, F120L/A266G, K438N/K500R, V79L/S509N, E175V/S257A/E368D/A469S, T71P/L374M, F280L/H440R, E89D/H440R, V58F/A328T/N532K, S119F/D342E/I449V, and K366N/G507S.
3. The isolated isoprene synthase variant of claim 1 or 2, wherein the variant has an N-terminal truncation.
4. The isoprene synthase variant of any one of the preceding claims, wherein the isoprene synthase variant is a variant of a plant isoprene synthase, selected from poplar, kudzu, aspen or English oak isoprene synthase.
5. The isolated isoprene synthase variant of any one of the preceding claims, wherein the variant comprises an amino acid sequence having at least 40% sequence identity to a wild-type isoprene synthase.
6. The isoprene synthase variant of claim 3, wherein the variant is selected from the group consisting of: an MEA variant (SEQ ID NO:122); an MSV variant (SEQ ID NO:124); an MVS variant (SEQ ID NO:126); an MTE variant (SEQ ID NO:128); an MNV variant (SEQ ID NO:130); a TRC (-3) variant (SEQ ID NO:136); a TRC (-4)

2009240505 02 Aug 2013

variant (SEQ ID NO:138); a TRC (-5) variant (SEQ ID NO:140); a TRC (-6) variant (SEQ ID NO:142); and a TRC (-7) variant (SEQ ID NO:144).

7. The isoprene synthase variant of claim 3, wherein the variant is a MET variant of *P. tremuloides* isoprene synthase (SEQ ID NO:146), or a MET variant of *P. trichocharpa* isoprene synthase (SEQ ID NO:148).
8. A host cell for producing isoprene comprising a heterologous polynucleotide sequence encoding an isoprene synthase variant of any one of claims 1 to 7 in operable combination with a promoter.
9. The host cell of claim 8, wherein the polynucleotide sequence is contained within a plasmid.
10. The host cell of claim 8, wherein the polynucleotide sequence is integrated into a chromosome of the host cell.
11. The host cell of any one of claims 8 to 10, wherein the host cell is selected from the group consisting of: gram-positive bacterial cells; gram-negative bacterial cells; filamentous fungal cells; yeast cells; *Escherichia sp.* (*E. coli*); *Pantoea sp.* (*P. citrea*); *Bacillus sp.* (*B. subtilis*); *Yarrowia sp.* (*Y. lipolytica*); and *Trichoderma* (*T. reesei*).
12. The host cell of any one of claims 8 to 11, wherein the host cell is cultured in a medium comprising a carbon source selected from the group consisting of glucose, glycerol, glycerine, dihydroxyacetone, yeast extract, biomass, molasses, sucrose, and oil.
13. The host cell of any one of claims 8 to 12, wherein the host cell further comprises a heterologous or native nucleic acid encoding an IDI polypeptide and/or a heterologous or native nucleic acid encoding a DXS polypeptide, optionally in combination with the native DXP pathway.
14. The host cell of any one of claims 8 to 12, wherein the host cell further comprises one or more nucleic acids encoding an IDI polypeptide and a DXS polypeptide.

2009240505 02 Aug 2013

15. The host cell of any one of claims 8 to 12, wherein the host cell comprises one vector encoding the isoprene synthase variant, an IDI polypeptide, and a DXS polypeptide.
16. The host cell of any one of claims 8 to 12, wherein the host cell further comprises a heterologous or native nucleic acid encoding a mevalonate (MVA) pathway polypeptide selected from the group consisting of an MVA pathway polypeptide from *Saccharomyces cerevisiae* and *Enterococcus faecalis*.
17. The host cell of any one of claims 8 to 12, wherein the host cell further comprises one or more nucleic acids encoding an MVA pathway polypeptide and a DXS polypeptide and wherein one vector encodes the isoprene synthase variant, the MVA pathway polypeptide, and the DXS polypeptide.
18. The host cell of claim 17, wherein the host cell further comprises one or more nucleic acids encoding a DXS polypeptide, an IDI polypeptide, or one or more DXP pathway polypeptides, and a MVA pathway polypeptide.
19. The host cell of any one of claims 8 to 12, wherein the host cell further comprises one or more nucleic acids encoding IDI and one or more MVA pathway polypeptide(s).
20. The host cell of claim 19, further comprising DXS.
21. A method of producing a host cell capable of producing isoprene, the method comprising introducing a heterologous polynucleotide sequence encoding an isoprene synthase variant of any one of claims 1 to 7 into the host cell.
22. A method of producing isoprene, comprising:
(a) culturing the host cells of any of claims 8 to 20 under suitable culture conditions for production of isoprene; and
(b) producing the isoprene.
23. The method of claim 22, further comprising:
(c) recovering the isoprene; and optionally
(d) polymerizing the isoprene.

02 Aug 2013

2009240505

24. The method of claim 22 or 23, wherein the host cell is selected from the group consisting of: gram-positive bacterial cells; gram-negative bacterial cells; filamentous fungal cells; yeast cells; *Escherichia sp.* (*E. coli*); *Pantoea sp.* (*P. citrea*); *Bacillus sp.* (*B. subtilis*); *Yarrowia sp.* (*Y. lipolytica*); *Trichoderma* (*T. reesei*); and *Saccharomyces* (*S. cerevisiae*).

25. The method of any one of claims 22 to 24, wherein the host cell is cultured in a medium that includes a carbon source selected from the group consisting of glucose, glycerol, glycerine, dihydroxyacetone, yeast extract, biomass, molasses, sucrose, and oil.

Figure. 1

Kudzu IspS DNA for *E. coli* (SEQ ID NO:1)

atgtgtgacgacctcttctcaatttactcagattaccgagcataattcccgctcgttccgcaaacct
gtggaatttcgaattcctgcaatccctggagaacgacctgaaagtggaaaagctggaggagaaaagcgaccaaacctgg
aggaagaagttcgcctgcatgatcaaccgtgtagacacccagccgctgtccctgctggagctgatcgacgatgtgcag
cgcttgggtctgacctacaaaatttgaaaaagacatcattaaaagccctggaaaacatcgtactgctggacgaaaacaa
aaagaacaaatctgacctgcacgcaaccgctctgtcttccgctcgtgctgcgtcagcacggttccgaggtttctcagg
atgtttttgagcgtttcaaggataaagaaggtgggttcagcgggtgaactgaaaggtgacgtccaaggcctgctgagc
ctgtatgaagcgtcttacctgggtttcgaggggtgagaacctgctggaggaggcgctaccttttccatcaccacac
gaagaacaacctgaaagaaggcattaataccaaggttgacagaacaagtgagccacgcctggaaactgccatatcacc
agcgtctgcacgctctggaggcacgttgggtccctggataaatacgaaccgaaaagaaccgcacaccagctgctgctg
gagctggcgaagctggattttaacatggtacagaccctgcaccagaaagagctgcaagatctgtcccgctgggtggac
cgagatgggcctggctagcaaaactggattttgtacgcgacccgctgatggaaagtattttctgggcaactgggtatgg
cgccagacccgcagtttgggtgaatgtcgcaaaagctgttactaaaatgtttgggtctgggtgacgatcatcgatgacgtg
tatgacgtttatggcactctggacgaactgcaactgttcaccgatgctgtagagcgctgggacgttaacgctattaa
caccctgccggactatatgaaactgtgtttcctggcactgtacaacaccggttaacgacacgctcctattctattctga
aagagaaaaggtcataacaacctgtcctatctgacgaaaagctggcggtgaactgtgcaaaagcctttctgcaagaggcg
aaatgggtccaacaacaaaattatcccggttttctccaagtacctggaaaacgccagcgtttcctcctccggtgtagc
gctgctggcgccgctcttacttttccgtatgccagcagcaggaagacatctccgaccacgcgctgcgttccctgaccg
acttccatgggtctgggtgcgttctagctgcgttatcttccgcctgtgcaacgatctggccacctctgcggcgagctg
gaacgtggcgagactaccaattctatcattagctacatgcacgaaaacgatggtagccagcgaggaaacaggcccgca
agaactgcgtaaactgatcgacgccgaatggaaaaagatgaatcgtgaacgcgttagcgactccaccctgctgccta
aagcgttcatggaaatcgcagttaacatggcacgtgttcccaactgcacctaccagtagggcgatgggtctgggtgcg
ccagactacgcgactgaaaaccgcacaaaactgctgctgattgacctttcccgattaacagctgatgtatgtcta
a

Figure 2

Kudzu IspS protein (SEQ ID NO:2)

MCATSSQFTQITEHNSRRSANYQPNLWNFEFLQSLENDLKVEKLEEKATKLEEEVRCMINRVDTQPLSLELIDDVQ
RLGLTYKFEKDIIKALENIVLLDENKKNKSDLHATALSFRLLRQHGFVSQDVFERFKDKEGGFSGELKGDVQGLLS
LYEASYLGFEGENLLEEARTFSITHLKNNLKEGINTKVAEQVSHALELPYHQRLHRLEARWFLDKYEPKEPHHQLLL
ELAKLDFNMVQTLHQKELQDLSRWWTMGLASKLDFVRDRLEMEVYFWALGMAPDPQFGECRKAVTKMFGLVTIIDDV
YDVYGTLDLQLFTDAVERWDVNAINLTPDYMKLCFLALYNTVNDTSYSILKEKGHNNLSYLTKSWRELCKAFLQEA
KWSNNKIIPAFSKYLENASVSSSGVALLAPSYFSVCQQQEDISDHALRSLTDFHGLVRSSCVIFRLCNDLATSAEL
ERGETTNSIISYMHENDGTSEEQAREELRKLIDAEWKMMNRERVSDSTLLPKAFMEIAVNMARVSHCTYQYGDGLGR
PDYATENRIKLLLDPPFINQLMYV

Figure 3

Poplar IspS DNA for *E. coli* (SEQ ID NO:6)

atgtgctctgtttctaccgagaacgtttccttcactgagacggaaaccgaggcacgtcgtagcgcgaactacgagcc
gaatagctgggactacgatttcctgctgtcttccgatactgacgaatctattgaggtgtacaaagacaaagcaaaga
aactggaggtgaagtgcgccgcgaaattaacaacgagaaagctgaattcctgactctgctggagctgatcgataac
gtacagcgctgggtctgggttaccgcttcgaatctgatatccgtcgcgcactggatcgtttcgttaagcagcggcgg
tttcgatggcgtgacaaaaacgagcctgcacgctaccgcgctgtccttcctgctcgtcgtcagcacggcttcgaag
tttctcaggaagcattctccggtttcaaagatcaaaacggtaacttcctggaaaacctgaaagaagacactaaggcg
atcctgagcctgtatgaggcaagctttctggccctggagggtgagaacatcctggatgaggcgcgctattcgccat
ctcccatctgaaagagctgtctgaagagaaaatcggttaaggaaactggcagagcagggttaatcacgcaactggaactgc
cgctgcctcgtcgtacccagcgtctggaggcgggttggtccatcgaagcgtaccgcaaaaaggaggatgctaaccag
gttctgctggaactggccatcctggactacaacatgatccagtcggtttaccagcgtgatctgctgaaacctcccg
ttgggtggcgcgctgtgggcctggcgacaaaactgcacttcgctaaggaccgcctgattgagtccttttactgggcag
tcggcgttgcggttcgaacctcagttctgactgccgtaacagcgttgcaaaaatgttcagcttcgttactattatc
gacgacatctacgacgtttacggctactctggacgagctggaactgtttaccgacgctgtcgaacgttgggatgttaa
cgccatcaacgatctgectgactacatgaaactgtgcttcctggcactgtataaacacgatcaacgaaattgcatacg
acaacctgaaagacaaagggtgaaaacatcctgccgtacctgactaaagcgtgggcggatctgtgtaacgcttttctg
caagaagcgaaatggctgtataacaaatccactccgaccttgacgattatttcggcaatgcctggaaatccagctc
tgcccgctgcaactgatcttcgcttattttgcggttggtccaaaacatcaaaaaggaggaaattgaaaacctgcaaa
aataccacgatatcattagccgtccttctcatatctttcgctgtgcaacgacctggcaagcgcgtccgcagagatc
gcacgtggcgaaaccgctaaactctgtttcctgctacatgcgcaccaagggcatttccgaagagctggcaaccgagag
cgtaatgaatctgatcgacgaaacctgtaagaaaatgaacaaagaaaaactgggtggctccctgttcgctaaaccgt
tcgtagagactgctattaacctggcagctcagagccactgcacctaccacaatgggtgacgcacatactagcccggt
gaactgactcgtaaacgtgtactgtctgttatcccgaaaccgattctgccgttcgaacgttaa

Figure 4

Poplar IspS protein (SEQ ID NO:7)

MCSVSTENVSFTEETETARRSANYEPNSWDYDFLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLLLELIDN
VQRLGLGYRFESDIRRALDRFVSSGGFDGVTKISLHATLSFRLLRQHGFEVSQEAFFSGFKDQNGNFLENLKEDTKA
ILSLYEASFLALEGENILDEARVFAISHLKELSEEKIGKELAEQVNHAEPLHRRRTQRLEAVWSIEAYRKKEDANQ
VLELAILDYNMIQSVYQRDLRETSRWRRVGLATKLHFAKDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTII
DDIYDVYGTLDLELELFTDAVERWDVNAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFL
QEAKWLYNKSTPTFDDYFGNAWKSSGPLQLIFAYFAVVQNIKKEEIENLQKYHDIISRPSHIFRLCNDLASASAEI
ARGETANSVSCYMRKGI SEELATESVMNLIDETCKMKNKEKLGGS LFAKPFVETAINLARQSHCTYHNGDAHTSPD
ELTRKRVLSVITEPILPFER

Figure 5


Kudzu Isoprene Synthase



Figure 6

Poplar Isoprene Synthase



**FIG. 7**




FIG. 8




FIG. 9

Figure 10A

plasmid MCM93 = pCR2.1-Kudzu (SEQ ID NO:22)

aagggcggaattctgcagatatccatcacactggcggcgctcgagcatgcatctagagggcccaattcgccctatag
tgagtcgtattacaattcactggccgctcggtttacaacgtcgtgactgggaaaaacctggcggtaccaccaactaatc
gccttgacgacacatcccccttgcgcagctggcgtaatagcgaagaggcccgaccgatcgcccttcccaacagttg
cgcagcctgaatggcgaatggacgcgcctgtagcggcgcattaaagcgcggcggtgtggtggttacgcgcagcgtg
accgctacacttgccagcgccttagcgcgcctccttgcgtttcttcccttcccttctcgccacgttcgcggcctt
tccccgtcaagctctaaatcgggggctccctttagggttccgatttagtgctttacggcacctcgaccccaaaaaac
ttgattagggtaggttcacgtagtgggccatcgccctgatagacgggttttgcgccttgacgttggagtccacg
ttctttaatagtgactcctgttccaaactggaacaacactcaaccctatctcggtctattcttttgattataagg
gattttgcccgtattcgccctattggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaaattc
agggcgcaagggctgctaaaggaagcggaacacgtagaaagccagtcgcgagaaacgggtgctgaccccggtgaatg
tcagctactgggtctatctggacaagggaaaaacgcaagcgcaaaagagaagcaggtagcttgagtgggcttacatgg
cgatagctagactggcggttttatggacagcaagcggaacgggaattgcccagctggggcgccctctggttaagggttg
gaagccctgcaaagtaaaactggatggcttcttgcgcgaaggatctgatggcgaggggatcaagatctgatcaag
agacaggatgaggatcggttcgcatgattgaacaagatggattgcacgcaggttctccggccgcttgggtggagagg
ctattcggtatgactgggcacacagacaatcggtgctctgatgccgcctgttccggtgtcagcgcagggcg
cccgttcttttgcgaagaccgacctgtccggtgcctgaatgaactgcaggacgagggcagcgcggctatcggtggc
tgccacgacggcggttccctgcgcagctgtgctcgacgttgcactgaagcgggaaggagctggctgctattgggc
gaagtgcggggcgaggtatcctgtcatccacacttgcctcctgcgagaaagtatccatcatggctgatgcaatgcg
gcggctgcatacgcttgatccggctacctgccattcgaccaccaagcgaaacatcgcatcgagcgagcacgtactc
ggatggaagccggctcttgctgatcaggatgatctggacgaagagcatcaggggctcgcccgacggcaactgttcgcc
aggctcaagggcgcatgcccgacggcgaggatctcgtcgtgacccatggcgatgcctgcttgcgcaatcatggt
ggaaaatggcgcttttctggattcatcgactgtggcgggtgggtgtggcggaaccgctatcaggacatagcgttgg
ctaccgctgatattgctgaagagcttggcgcggaatgggtgacccgcttctcgtgcttaccgtatcgccgctccc
gattcgacgcgcacgccttctatcgcccttcttgacgagttcttctgaattgaaaaaggaagatgagttatgaa
catttccgtgtcgcccttattcccttttgcggcatttgccttctgttttgcctcaccagaaacgctggtgaa
agttaaagatgctgaagatcagttgggtgcacgagtggttacatcgaaactggatctcaacagcggttaagatccttg
agagtttgcgccgaagaacggtttccaatgatgagcactttaaaagtctcgtatgtggcgcggtattatcccgt
attgacgcgggcaagagcaactcggtcgccgcatacactattctcagaatgacttgggtgagtagtaccagtcac
agaaaagcatcttacggatggcatgacagtgaagaaattatgcagtgctgccataacctagtgataaactgcgg
ccaacttacttctgacaacgatcgaggaccgaaggagctaaacgctttttgcacaacatgggggatcatgtaact
cgcccttgatcggttgggaaccggagctgaatgaagccataccaaacgacgagcgtgacaccacgagcctgtagcaat
ggcaacaacgcttgcgcaactattaactggcgaactacttactctagcttcccggaacaattaatagactggatgg
aggcgataaaagtgcaggaccacttctgcgtcggcccttccggctggctgggtttatgctgataaatctggagcc
ggtgagcgtgggtctcgcggtatcattgcagcactggggccagatggttaagccctcccgtatcgtagttatctacac
gacggggagttaggcaactatggatgaacgaatatagacagatcgctgagataggcgctcactgattaaagcatggg
aactgtcagaccaagttactcatatatactttagattgatttaaaacttcaatttttaatttaaaaggatctagggtg
aagatcctttttgataatctcatgacaaaaatcccttaacgtgagtttctggttccactgagcgtcagaccccgtaga
aaagatcaaaggatcttcttgagatccttttttctgcgcgtaattctgctgcttgcaaacaaaaaaaccaccgctac
cagcgggtgggtttgttgcgggatcaagagctaccaactcttttccgaaggtaactggcttcagcagagcgagata
ccaaataactgttcttctagttagcgttagttaggccaccacttcaagaactctgtagcaccgctacatccctcgc
tctgctaatcctgttaccagtggtgctgccagtgccgataagtcgtcttaccgggttgagctcaagacgatagt
taccggataaaggcgacgggtcggtgctgaacggggggttcgtgcacacagccagcttggagcgaacagctacacc
gaactgagataacctacagcgtgagctatgagaaagcgccacgcttccgaaggagaaaagcgagcaggtatccggt
aagcggcaggggtcggaacaggagagcgcacagggagcttccaggggaaaacgcctggtatctttatagtcctgtcg
ggtttcgccacctctgacttgagcgtcgattttgtgatgctcgtcagggggcgagcctatggaaaaacgcccagc
aacggcgcccttttacgggtcctggccttttgcgtgcttgcacatgttcttccctgcgttatccctgattc
tgtggataaacgctattaccgctttgagtgagctgataccgctcgccgcagccgaacgaccgagcgcagcagtgag
tgagcgaggaagcggaagagcgcccaatacgcaaaccgctctcccgcgcggttggccgattcattaatgcagctgg
cacgacaggtttcccgactggaaaagcgggcagtgagcgcaacgcaattaatgtgagttagctcactcacttaggcacc
ccaggtttacactttatgcttccggctcgtatgttggtggaattgtgagcggataaacaattcacacaggaaca
gctatgacctgattacgcaagcttggtagcagctcggtatccnctagtaacggccgagcgtgtgctggaattcgc
ccttgatcatgcatcgccttaggaggtaaaaaaacatgtgtgcgaccttctcacttactcagattaccgagc
ataattcccgctgcttccgcaaacatcagccaaacctgtggaatttcaattcctgcaatccctggagaacgacctg

Figure 10B

aaagtggaaaagctggaggagaaagcgaccaaactggaggaagaagttcgctgcatgatcaaccgtgtagacacca
gccgctgtccctgctggagctgatcgacgatgtgcagcgctgggtctgacctacaaatttgaaaaagacatcatta
aagccctggaaaacatcgtaactgctggacgaaaacaaaaagaacaaatctgacctgcacgcaaccgctctgtcttc
cgtctgctgctgcgcagcacggtttcgagggtttctcaggatgtttttgagcgtttcaaggataaagaaggtggtttcag
cggtgaaactgaaaggtgacgtccaaggcctgctgagcctgtatgaagcgtcttacctgggtttcgaggggtgagaacc
tgctggaggaggcgcgtaacctttccatcacccacctgaagaacaacctgaaagaaggcattaataccaaggttgca
gaacaagtgagccacgccctggaactgccatataccacagcgtctgcaccgtctggaggcacgttggttccctggataa
atacgaaccgaaagaaccgcatcaccagctgctgctggagctggcgaaagctggattttaacatggtacagaccctgc
accagaaaagagctgcaagatctgtcccgtgggtggaccgagatgggcctggctagcaaaactggattttgtacgcgac
cgctgatggaagtttatcttctgggcactgggtatggcgccagaccgcagtttggtgaatgtcgcaaagctgttac
taaaatgtttggtctggtgacgatcatcgatgacgtgtatgacgtttatggcactctggacgaactgcaactggtca
ccgatgctgtagagcgctgggacgttaacgctattaacaccctgccggactatatgaaactgtgtttcctggcactg
tacaacaccgttaacgacacgtcctattctattctgaaagagaaaggtcataacaacctgtcctatctgacgaaaag
ctggcgtgaactgtgcaaagcctttctgcaagaggcgaaatggtccaacaacaaaattatcccggctttctccaagt
acctggaaaacgcccagcgtttcctcctccggtgtagcgtctggcgccgtcttactttccgtatgccagcagcag
gaagacatctccgaccacgcgctgcgttccctgaccgacttccatggtctggtgcgttctagctgcgttatcttccg
cctgtgcaacgatctggccacctctgcggcgagctggaacgtggcgagactaccaattctatcattagctacatgc
acgaaaaagatggtaccagcgaggaaacaggcccgcgaagaactgcgtaaactgatcgacgccgaatggaaaaagatg
aatcgtgaacgcgttagcgactccacctgctgcctaaagcgttcatggaaatcgagttaacatggcacgtgtttc
ccactgcacctaccagtatggcgatggtctgggtcgcccagactacgcgactgaaaaccgcatcaaactgctgctga
ttgaccctttcccgatataccagctgatgtatgtctaactgcagggatccgtcgaccg


**FIG. 11**

Figure 12A

pET24D-Kudzu (SEQ ID NO:23) Kudzu IspS ORF 48-1742 (complementary)
gtgcgccgcaagcttgcgacggagctcgaattcggatccctgcagttagacatacatcagctgggttaatcgggaa
agggccaatcagcagcagcttgcggttttcagtcgcgtagctctggcgacccagaccatcgccatactggtagg
tgcagtgaggaaacacgtgccatgttaactgcgatttccatgaacgctttaggcagcagggtaggagtcgctaacgcgt
tcacgattcatcttttccattcggcgctcgatcagtttacgcagttcttcgcgggctgttccctcgctggtagccatc
gttttcgtgcatgtagctaatgatagaattggtagtctcgccacgttccagctccgcccagagggggccagatcgt
tgcacaggcggaagataaacgcagctagaacgcaccagaccatggaagtgcgtcagggaacgcagcgctggtcggag
atgtcttccctgctgctggcatacggaagtaagacggcgccagcagcgtacaccggaggaggaaacgcgtggcggt
ttccaggtacttggagaaagccgggataattttgttggtagcatttcgcctcttcgcagaaaggctttgcacagtt
cacgccagcttttcgtcagataggacaggttggtatgacctttctcttcagaatagaataggacgtgtcgttaacg
gtgtgtacagtgccaggaaacacagtttcatatagtcggcaggggtgttaatagcgttaacgtcccagcgctctac
agcatcggtgaacagttgcagttcgtccagagtgccataaacgtcataacagtcacatcgatgatcgtcaccagacca
acatttttagtaacagctttgcgacattcaccaaactcggggtctggcgccataccagtgcccagaaataaacttcc
atcaggcggtcgcgtagaaaaatccagtttgctagccaggcccatctcggtccaccagcgggacagatcttcgagctc
tttctggtagcgggtctgtaccatgttaaaatccagcttcgcagctccagcagcagctggtagcgggtcttctcg
gttcgtattttatccaggaaaccaacgtgcctccagacgggtgcagacgctggtagatggcagttccaggcggtggctc
acttgttctgcaaccttgggtatataatgccttcttccaggttggtcttcaggtgggtgatggaaaaggtagcgcctc
ctccagcaggttctcaccctcgaaacccaggttaagacgcttcatacaggctcagcaggccttgagcgtcaccttca
gttcacgcgtgaaacacaccttcttctccttgaaacgctcaaaaaacatcctgagaaacctcgaaacacgtgctgacgc
gcagacggaaaagacagagcggttgcgtgcaggtcagatttggtcttttcttcgtccagcagtagcgttttcc
cagggtcttaatgatgtcttttcaaatttgtaggtcagacccaggcgctgcacatcgtcgatcagctccagcaggg
acagcggtgggtgtctacacggttgatcatgcagcgaacttcttccctcagtttggtcgcttctcctccagcttt
tccactttcaggtcgttctccagggttgaggaattcgaaattccacaggttggtgatagttgcggaacgacg
ggaattatgctcggtaatctgagtaaatgagaagaggtcgacacatggtatatctccttctaaagttaaacaaa
attattctagaggggaattgttatccgctcacaattcccctatagtgagtcgtattaatttcgcgggatcgagatc
tcgatcctctacgcgcgagcgcgttcggtcgccatcaccggcgccacaggtgcgggtgctggcgcttatcgcgga
catcaccgatggggaagatcgggctcgccacttcgggtcatgagcgcttggttcggcggtgggtatggtagcggcc
cggtggccgggggactgtggcgccatctccttgcatgcacattccttcggcgcggtgctcaacggcctcaac
ctactactgggtgcttccctaatgcaggagtcgcataaggagagcgtcgagatcccggacaccatcgaatggcgca
aaaccttctcgcggtatggcatgatagcggccggaagagagtcgaattcagggtgggtgaatgtgaaaccagtaacgtta
tacgatgtcgagagatagcgggtgtctcttatcagaccgtttcccgcggtggtagaacaggccagccacgtttctgcg
gaaaacgcgggaaagtggaagcggcgatggcgagctgaattacattcccaaccggtggcacaacacacgtggcg
gcaaacagtcgttgctgattggcggtgccacctccagctcgccctgcacgcgcgtcgcaaatgtcgcgcggtat
aatctcgcgccgatcaactgggtgccagcgtgggtgtcgatggtagaacgaagcggcgctcgaagcctgtaaagc
ggcggtgcacaatcttctcgcgcaacgcgtcagtggggtgatcattaaactatccgctggatgaccaggatgccattg
ctgtggaagctgcctgcactaatgttccggcggtatttcttgatgtctctgaccagacaccatcaacagttattt
ttctcccatgaagacggtacgcgactggcggtggagcatctggtagcattgggtcaccagcaaatcgcgctgttagc
gggcccattaaagtctgtctcgcgcgctctgcgtctggctggctggcataaatatctcactcgcaatcaaatcagc
cgatagcggaacgggaaggcgactggagtccatgtccggttttcaacaaaccatgcaaatgctgaatgagggcatc
gttccactgcgtagctgggtgccaaacgatcagatggcgctggcgcaatgcgcgccattaccaggtccgggtgcg
cggtgggtgcggtatctcggtagtggtatcagacgataccgaagacagctcatgttatatcccgccgttaaccacca
tcaaacaggatttttcgctgctggggcaaacacgcgtggaaccgcttgctgcaactctctcaggggccaggcggtgaag
ggcaatcagctgttgcccgctctcactggtgaaaagaaaaaccacctggcgcccaatcgcgaacacgcctctccccg
cgcttgccgatctcattaatgcagctggcagcagaggtttcccgactggaaagcgggcagtgagcgcaacgcaatt
aatgtaagttagctcactcattaggcaccgggtctcgaccgatgcccttgagagccttcaaccagtcagctcctt
ccggtggcgcggggcatgactatcgtcgcgcaacttatgactgtcttctttatcatgcaactcgtaggacaggtgc
cggcagcgctctgggtcattttcggcgaggaccgctttcgctggagcgcgacgatgatcggtgtcgcttgcggtta
ttcggaattctgcacgcctcgctcaagccttcgtcactgggtcccgcaccaaacgtttcgcgagaaagcaggccat
tatcgccggcatggcgggccacgggtgcgcgatgatcgtgctcctgtcgttgaggacccggctaggctggcggggtt
gccttactgggttagcagaatgaatcaccgatcgcgagcgaacgtgaagcgactgctgctgcaaacgctctgcgacc
tgagcaacaacatgaatggcttcggtttccgtgtttcgttaaagctcggaacgcgggaagtcagcgccctgcaccat
tatgttcggatctgcacgcaggtgctgctggctaccctgtggaacacctacatctgtattaacgaagcgctggc
attgacctgagtgatttttctgtgtcccgccgcacccataccgccagttgtttacctcacaacgttccagttaac
cgggcatgttcatcatcagtaaacccgtatcgtgagcatcctctcgtttcatcggtatcattaccccatgaacag

Figure 12B

aaatcccccttacacggagggcatcagtgaccaaacaggaaaaaacgccttaacatggcccgctttatcagaagcc
agacattaacgcttctggagaaactcaacgagctggacgcggatgaacaggcagacatctgtgaatcgcttcacgac
cacgctgatgagctttaccgcagctgcctcgccgctttcggtgatgacggtgaaaaacctctgacacatgcagctccc
ggagacgggtcacagcttgtctgtaagcggatgcccgggagcagacaagcccgctcagggcgcgctcagcgggtgttggcg
ggtgtcggggcgagccatgacccagtcacgtagcgatagcggagtgatactggcttaactatgcggcatcagagc
agattgtactgagagtgacccatataatgcggtgtgaaataaccgcacagatgcgtaaggagaaaaataccgcatcaggc
gctcttcgcttctcgcctcactgactcgctgcgctcggtcgttcggctcgccgagcggtatcagctcactcaaag
gcggtaataggttatccacagaatcaggggataacgcaggaaagaacatgtgagcaaaaggccagcaaaaggccag
gaaccgtaaaaaggccgcttgctggcgttttccataggctccgccccctgacgagcatcacaaaaatcgacgct
caagtgcagaggtggcgaaacccgcagaggactataagataaccaggcggttccccctggaagctccctcgtgcgctct
cctgttcgacccctgcgcttaccggatacctgtccgcttctcccttcgggaagcggtggcgctttctcatagctc
acgctgtaggtatctcagttcggtgtaggtcgttcgctccaagctgggctgtgtgcagaaacccccggttcagccc
accgctgcgcttatccggttaactatcgtcttgagtcacacccggttaagacacgacttatcgccactggcagcagcc
actggtaacaggattagcagagcgaggtatgtaggcggtgtacagagttcttgaagtgggtggcctaactacggcta
cactagaaggacagttatgtgtatctgcgctcgtcgtgaagccagttaccttcggaaaaagagttggtagctcttgat
ccggcaaaacaaaccacgctggtagcggtgggtttttgtttgcaagcagcagattacgcgcagaaaaaaaggtatct
caagaagatcctttgatctttctacggggtctgacgctcagtggaacgaaaactcacgttaagggattttggtcat
gaacaataaaaactgtctgcttacataaaacagtaatacaaggggtgttatgagccatattcaacgggaaacgtcttg
tctaggccgcgattaaattccaacatggatgctgatttatatgggtataaatgggctcgcgataatgtcgggcaatc
agggtgcgacaatctatcgattgtatgggaagcccgatgcgccagagttgttctgaaacatggcaaaaggtagcggtg
ccaatgatgttacagatgagatgggtcagactaaactggctgacggaatttatgcctcttcggacctcaagcatttt
atccgtactcctgatgatgcattggttactcaccactgcgctcccgggaaaaacagcattccagggtattagaagaata
tcttgattcaggtgaaaaatattgttgatgcgctggcagtggtctcgcgcgggttgcatctgattcctgtttgtaat
gtccttttaacagcgatcgctatttcgtctcgtcagcgcaatcacgaatgaataacgggtttggttgatgcgaggt
gattttgatgacgagcgtaaatggctggcctgttgaacaagcttggaagaaatgcataaaacttttgccattctcacc
ggattcagtcgctcactcatggtgattttcacttgataaccttatttttgacgaggggaaataataggttgattg
atgttgagcagagtcggaatcgacagccgataccaggatcttgccatcctatggaactgcctcggtgagttttctcct
tcattacagaaacggtttttcaaaaatattggtattgataatcctgatagaataaattgcagtttcatattgatgct
cgatgagttttctaaagaattaattcatgagcggatataattgaatgtatttagaaaaataaacaataaggggtt
ccgcgcacatttccccgaaaagtgccacctgaaattgtaaacgttaataattttgttaaaattcgcgtaaattttg
ttaaatcagctcatttttaaccaataggccgaaatcgcaaaaatcccttataaatcaaaagaatagaccgagatag
gggtgagtggttccagtttgaacaagagtcactattaaagaacgtggactccaacgtcaaaggcgaaaaaac
gtctatcagggcgatggccactacgtgaaccatcacctaatcaagtttttggggtcgaggtgccgtaaaagcact
aatcggaacctaaaaggagccccgatttagagcttgacggggaaagccggcgaaacgtggcgagaaaggaaaggga
agaaagcgaaaggagcgggcgctagggcgctggcaagttagcgggtcacgctgcgcgtaaccaccacacccgcgcg
cttaatgcgcgctacagggcgcgctccattcgccaatccggatatagtctcctttcagcaaaaaaccctcaag
accggttagaggcccaagggttatgctagtatttgctcagcggtggcagcagccaactcagcttctttcgggc
tttgtagcagcggtatcagtggtggtggtggtggtgctcga




FIG. 13




FIG. 14

Figure 15




Figure 16




Figure 17

Amino acid sequence of 6XHis N-terminally tagged *P. alba* IspS (SEQ ID NO:118)

MRGSHHHHHHGMASMTGGQMGRDLYDDDDKDHPFTMRCSVSTENVSFTETETEARRSANYEPN
SWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLELIDNVQRLGLGYRFESDI
RGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFVSQEAFFSGFKDQNGNFLENLKEDIKAI
LSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVNHAELEPLHRRTORLEAVWS
IEAYRKEDANQVLELAILDYNMIQSVYQRDLRETSRWWRRVGLATKLHFARDRLIESFYWAV
GVAFEPQYSDCRNSVAKMFSFVTIIDDIYDVYGTLDLELELFTDAVERWDVNAINDLPDYMKLCF
LALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNKSTPTFDDYFGNAWKSSS
GPLQLVFAYFAVVQNIKKEEIEENLQKYHDTISRPSHIFRLCNDLASASAEIARGETANSVSCYM
RTKGISEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLARQSHCTYHNGDAHTSPD
ELTRKRVLSVITEPILPFER

Figure 18A

Nucleotide sequence of pDu27 (SEQ ID NO:119)

aaggcgagctcaacgatccggctgctaacaaagcccgaaaggaagctgagttggctgctgccac
ccgtgagcaataactagcataaacccttggggcctctaaacgggtcttgaggagtttttggct
gaaaggaggaactatatccggatatcccgcaagaggcccggcagtagccggcataaccaagccta
tgcttacagcatccagggtgacgggtgccgaggatgacgatgagcgcatgttagatttcataca
cgggtgcctgactgcttagcaatttaactgtgataaactaccgcattaaagcttatcgatgata
agctgtcaaacatgagaattaattcttgaagacgaaaggcctcgtgatacgcctattttata
ggttaatgtcatgataataatggtttcttagacgtcaggtggcacttttcggggaaatgtgcgc
ggaaccctatttgtttatttttctaaatacattcaaataatgtatccgtcatgagacaataac
cctgataaatgcttcaataatattgaaaaaggaagagtattgattgaacaagatggattgcacgc
aggttctccggccgcttgggtggagaggctattcggctatgactgggcacaactgacaatcggc
tgcctgatgccgcgctgttccggctgtcagcgcaggggcgcccggttcttttgcgaagaccg
acctgtccgggtgccctgaatgaactgcaggacgaggcagcgcggtatcgtggctggccacgac
ggcggttcttgcgcagctgtgctcgacgttgtcactgaagcgggaaggactggctgctattg
ggcgaagtgcgggggcaggatctcctgtcatctcaccttgcctcctgccgagaaagtatccatca
tggctgatgcaatgcggcggtgcatacgccttgatccggctacctgccattcgaccaccaagc
gaaacatcgcatcgagcgggcacgtactcggatggaagccgggtcttgcgatcaggatgatctg
gacgaagagcatcaggggtcgcgccagccgaactgttcgccaggctcaaggcgcgatgcccg
acggcgaggatctcgtcgtgacacatggcgatgcctgcttgcgaatatcatgggtggaaatgg
ccgcttttctggattcatcgactgtggcgggtgggtgtggcggaaccgctatcaggacatagcg
ttggctaccgctgatattgctgaagagcttggcggcgaatgggctgaccgcttctcgtgcttt
acggtatcgccgctcccgattcgcagcgcacgccttctatcgcttcttgacgagttcttctg
agcgggactctggggttcgaaatgaccgaccaagcgacgcctaaactgtcagaccaagttactc
atatatacttttagattgattttaaacttcatttttaatttaaaggatctaggtgaagatcctt
tttgataatctcatgacaaaaatcccttaacgtgagttttcgttccactgagcgtcagaccccg
tagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaaac
aaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttccg
aaggttaactggcttcagcagagcgcagataccaaataactgtccttctagttagccgtagttag
gccaccacttcaagaactctgtagcaccgcctacatacctcgctctgctaatcctgttaccagt
ggctgctgccagtggcgataagtcgtgtcttaccgggttggactcaagacgatagttaccggat
aaggcgcagcggtcgggctgaacggggggttcgtgcacacagcccagcttggagcgaacgacct
acaccgaactgagatacctacagcgtgagctatgagaaaagcggccagcttcccgaaggagaaa
ggcggacaggtatccggtgaacgggcagggtcggaacaggagagcgcacgagggagcttccaggg
ggaaacgcctggtatctttatagtctgtcggtttcgccacctctgacttgagcgtcgatttt
tgtgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgggt
cctggccttttgcgtggccttttgcctcacatgttctttcctgcgttatccctgattctgtggat
aaccgtattaccgcctttgagtgcgtgataccgctcgccgcagccgaacgaccgagcgcagcg
agtgcgtgagcgcaggaagcgggaagagcgcctgatgcgggtattttctccttacgcattctgtgcgg
tatttcacaccgcaatggtgcactctcagtacaatctgctctgatgccgcatagttaagccagt
atacactccgctatcgctacgtgactgggtcatggctgcgccccgacaccgcgaacaccgcgt
gacgcgccctgacgggcttgtctgctcccggcatccgcttacagacaagctgtgaccgtctccg
ggagctgcatgtgtcagaggttttaccgctcatcaccgaaacgcgcgagggcagctgcggtaag
ctcatcagcgtggctcgtgaagcattcacagatgtctgcctgttcatccgcgtccagctcgttg
agtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaagggcggttttt

Figure 18B

Nucleotide sequence of pDu27

cctgtttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgataccgat
gaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactggaacgt
tgtgagggtaaacaactggcggtatggatgcggcgggaccagagaaaaatcactcagggtcaat
gccagcgcttcgttaatacagatgtaggtgtccacagggtagccagcagcatcctgcgatgca
gatccggaacataaatggtgcagggcgctgacttccgcgtttccagactttacgaaacacggaaa
ccgaagaccattcatgttggtgctcaggtcgcagacgttttgcagcagcagtcgcttcacgttc
gctcgcgtatcgggtgattcattctgctaaccagtaaggcaaccccgcagcctagccgggtcct
caacgacaggagcacgatcatgcgcacccgtggccaggacccaacgctgcccagatgcgcgc
gtgcggtgctggagatggcgacgcgatggatgttctgccagggttggtttgcgcattca
cagttctccgcaagaattgattggctccaattcttggagtgggaatccgttagcgaggtgccg
ccggttccattcaggtcgaggtggcccggtccatgcaccgcgacgcaacgcggggaggcaga
caaggtatagggcggcgcctacaatccatgccaacccgttccatgtgctcgccgagggcgcata
aatcgccgtgacgatcagcgggtccaatgatcgaagttagggtggttaagagccgcgagcgatcct
tgaagctgtccctgatggtcgtcatctacctgcctggacagcatggcctgcaacgcgggcaccc
cgatgccgcgggaagcgagaagaatcataatggggaaggccatccagcctcgcgctcgcaacgc
cagcaagacgtagcccagcgctcgcccgccatgccggcgataatggcctgcttctcgccgaaa
cgtttggtggcgggaccagtgcgaaggttgagcgagggcggtgcaagattccgaataccgcaa
gcgacaggccgatcatcgtcgcgctccagcgaaaagcggtcctcgccgaaaatgacccagagcgc
tgccggcacctgtcctacgagttgcatgataaagaagacagtcataagtgcggcgacgatagtc
atgccccgcgcccaccggaaggagctgactgggttgaaggctctcaagggtcaggtcgagatc
ccggtgcctaatgagtgcgttaacttacattaattgcgttgcgctcactgcccgtttccagtc
gggaaacctgtcgtgccagctgcattaatgaatcgcccaacgcgcggggagaggcggtttgcgt
attgggcgccagggtggtttttctttcaccagtgcgacgggcaacagctgattgcccttcacc
gcctggccctgagagagttgcagcaagcggtccacgctgggtttgcccagcaggcgaaaatcct
gtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcgtcgtatcccactac
cgagatatccgcaccaacgcgcagcccggactcggtaatggcgcgcatcgcgccagcgccatc
tgatcgttggcaaccagcatcgagtggaacgatgccctcattcagcatttgcatggtttgtt
gaaaaccggacatggcactccagtcgccttccggttccgctatcggtgaatttgattgcgagt
gagatatttatgccagccagccagacgcgacgcgcgagacagaacttaatgggcccgttaac
agcgcgatttgctggtgacccaatgcgaccagatgctccacgcccagtcgcgtaccgtcttcat
gggagaaaaataatactgttgatgggtgtctgggtcagagacatcaagaaataacgcgggaacatt
agtgcaggcagcttccacagcaatggcatcctgggtcatccagcggatagttaatgatcagccca
ctgacgcgttgcgcgagaagattgtgcaccgcgcgtttacaggcttcgacgcgcgttcgttcta
ccatcgacaccaccacgctggcaccagttgatcggcgcgagatttaatcgccgcgacaatttg
cgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacgactgtttgcccgcc
agttgttgcacacgcggttgggaatgtaattcagctccgccatcgccgcttccacttttccc
gcgttttcgcagaaacgtggctggcctggttcaccacgcgggaaacgggtctgataagagacacc
ggcatactctgcgacatcgtataacgttactggtttcacattcaccacctgaattgactctct
tccgggcgctatcatgccataccgcgaaagggttttgccgcatcogatggtgtccgggatctcga
cgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggttgaggccgttgagca
ccgccgcgcaaggaatggtgcattgcaaggagatggcgcccaacagtcccccggccacggggcc
tgccaccataccacgcgcaaaacaagcgctcatgagcccgaagtggcgagcccgatcttcccca

Figure 18C

Nucleotide sequence of pDu27

tccggtgatgtcggcgatataaggcgccagcaaccgcacctgtggcgccggtgatgccggccacga
tgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatacgactcactataggggaa
ttgtgagcggataacaattcccctctagaaataattttgtttaactttaagaaggagatataca
tatgcgggggttctcatcatcatcatcatcatggtagctgactggtggacagcaaatg
ggtcggggtctgtacgacgatgacgataaggatcatcccttcaccatgcgttgtagcgtgtcca
ccgaaaatgtgtctttcaccgaaactgaaaccgaagctcgtcgttctgcgaactacgaacctaa
cagctgggactatgattacctgctgtcctccgacacggacgagtcctcgaagtatacaaagac
aaagcgaaaaagctggaagccgaagttcgtcgcgagattaataacgaaaaagcagaatttctga
ccctgctggaactgattgacaacgtccagcgctgggctgggttaccgtttcgagctctgat
ccgtggtgcgtggatcgcttcgtttcctccggcggttcgatgcggtaaccaagacttccctg
cacggtacggcactgtctttccgtctgctgcgtcaacacggttttgaggtttctcaggaagcgt
tcagcggcttcaaagacccaaaacggcaacttccctggagaacctgaaggaagatatcaaagctat
cctgagcctgtacgaggccagcttccctggctctggaaggcgaaaacatcctggacgagcggaag
gttttcgcaatctctcatctgaaagaactgtctgaagaaaagatcggtaaagagctggcagaac
aggtgaaccatgcactggaactgccactgcacgcgtaactcagcgtctggaagcagtatggtc
tatcagggcctaccgtaaaaaggaggacgcgaatcaggttctgctggagctggcaattctggat
tacaacatgatccagtctgtataccagcgtgatctgcgtgaaacgtcccgttgggtggcgtcgtg
tgggtctggcgacccaaactgcactttgctcgtgaccgcctgattgagagcttctactgggcgt
gggtgtagcattcgaaccgcaatactccgactgccgttaactccgtcgcaaaaatgttttcttc
gtaaccattatcgacgatatactacgatgtatacggcaccctggacgaactggagctgtttactg
atgcagttgagcgttgggacgtaaacgccatcaacgacctgccggttacatgaaactgtgctt
tctggctctgtataaactattaacgaaatcgctacgacaacctgaaagataaagggtgagaac
atcctgccgtatctgaccaaagcctgggtgacctgtgcaacgctttcctgcaagaagccaagt
ggctgtacaacaaatctactccgaccttgacgactacttcggcaacgcatggaaatcctcttc
tggcccgctgcaactgggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatc
gaaaacctgcaaaaataccatgacaccatctctcgtccttcccatatcttccgtctgtgcaatg
acctggctagcgcgtctgcggaaattgcgcgtggtgaaaccgcaaatagcgtttcttggtacat
gcgcactaaaggatctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacc
tgaaaaaagatgaacaaggaaaaactgggtggtagcctgttcgcgaaaccgttcgtggaaaccg
cgatcaacctggcacgtcaatctcactgcacttatcataacggcgacgcgcataacctctccgga
tgagctgaccgcgaaacgcgttctgtctgtaatcactgaaccgattctgccgtttgaacgctaa

Figure 19

Amino acid sequence of full length *P. alba* IspS in *P. alba* (SEQ ID NO:120)

MRCSVTENSFTETETEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINN
EKAEFLTLLELIDNVQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGF
EVSQEAFSGFKDQNGNFLENLKEDIKAILSLYEASFLAEGENILDEAKVFAISHLKELSEEKI
GKELAEQVNHAELPLHRRTQRLEAVSIEAYRKKEDANQVLLELAILDYNMIQSVYQRDLRET
SRWWRRVGLATKLHFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDIYDVYGTLD
ELELFTDAVERWDVNAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNA
FLQEAKWLYNKSTPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEEIENLQKYHDTISRPSH
IFRLCNDLASAEIARGETANSVSCYMRTKGISELATESVMNLIDETWKKMNKEKLGGSLFA
KPFVETAINLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 20A

Nucleotide sequence of *P. alba* pET24a (SEQ ID NO:121)

tggcgaatgggacgcgcctgtagcggcgcattaaagcgcggcgggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccttagcgcgcctcctttcgctttcttcccttctttctcg
ccacgttcgcggctttccccgtcaagctctaaatcgggggctccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaaacttgattagggtgatggtcacgtagtggccatcg
ccctgatagacggtttttcgcccttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgattataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcagggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaatatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaataagg
tatcaagtgagaaatcaccatgagtgcgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatcaacc
aaaccgttattcattcgtgattgcgcctgagcgcgagcaaatacgcgatcgctgttaaaaggac
aattacaaacaggaatcgaatgcaaccggcgcaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaataacctggaatgctgttttcccggggatcgcagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcatggcaacgctacctttgccatgtttcag
aaacaactctggcgcacatgggcttcccatacaatcgatagattgtcgacactgattgcccagaca
ttatcgcgagccatttatacccatataaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttcccggtgaatatggctcataaacacccttgtattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcgttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggtaactggcttcagcagagcgcagataccaaataactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgctacatacctcgctctgctaatacctgttacca
gtggctgctgccagtggcgataagtctgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgcagcggctcgggctgaacggggggttcgtgcacacagcccagcttggaagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaaagcgccacgcttcccgaagggaga
aaggcggacaggtatccggtgaagcggcaggggtcggaacaggagagcgcacaggggagcttccag
ggggaacgcctggtatctttatagtctgtcgggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcaggggggaggagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatcccctgattctgtgg
ataaccgtattaccgcctttgagtgcgtgataccgctcgccgcagccgaacgaccgagcgcag
cgagtcagtgcgaggaagcgaagagcgcctgatgcggtattttctccttacgcatctgtgc
ggtatttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacaccgcgaacac
ccgctgacgcgcctgacgggttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cggttaggtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaaggggcgt
tttttctgtttggtcactgatgcctccgtgaaggggatttctgttcatgggggtaataata

Figure 20B

Nucleotide sequence of P. alba pET24a

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacactggcggatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggttccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcaggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcaggtgcagacgttttgcagcagcagtcgcttca
cggttcgctcgcgatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgatacgcgcacccgtggggccgcatgccggcgataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtgacgaaggcttgagcgagggcggtgcaagat
tccgaataccgcaagcgacaggccgatacgtcgcgctccagcgaaagcggtcctcgccgaaa
atgacccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagt
cgggcagcatagtcagtcggcgcccgccacgggaaggagctgactgggttgaaggctctcaagg
catcggtcgagatcccggtgctaatgagtgagctaacttacattaattgcgttgcgctcactg
ccgctttccagtcgggaaacctgtcggtccagctgcattaatgaatcgccaacgcgcgggga
gaggcggtttgcgtattggcgccagggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttcaccgctggcctgagagagttgcagcaagcggtccacgctgggttggccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtctcggtatcg
tcgtatccactaccgagatatcgccaccaacgcgcagcccgactcggtaatggcgcgcattg
cgccagcgccatctgatcggtggcaaccagcatcgagtggaacgatgccctcattcagcat
ttgatggtttgttgaaaaccggacatggcactccagtcgccttccggttccgctatcggtga
atattgatcgagtgagatatttatgccagccagccagacgcagacgcgcgagacagaaactta
atgggcccgcataacagcgcatattgctggtgacccaatgcgaccagatgctccacgccagtcg
cgtaccgtcttcatgggagaaaataactgttgatgggtgtctggtcagagacatcaagaaat
aacgcgggaacattagtgaggcagcttccacagcaatggcatcctggtcatccagcggatagt
taatgatcagccactgacgcgttgccgcgagaagattgtgcaccgcccgtttacaggcttcgac
gcccgttctgttaccatcgacaccaccagctggcaccagttgatcggcgcgagatttaac
gccgcgacaatttgcgacggcgctgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgcacgcgggttgggaatgtaattcagctccgccatcgccgc
ttccactttttcccgcttttcgcagaaacgtggtggcctggttaccacgcgggaaacgggtc
tgataagagacaccggcatactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggcgttgagcaccgcccgcgcaagggaatggtgcatgcaaggagatggcgcccaacagtc
ccggccacggggcctgccaccataccacgcggaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcgcgatataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatacgac
tcactataggggaattgtgagcggataacaattccctctagaaataattttgtttaactttaa
gaaggagatatacatatgcgttgtagcgtgtccaccgaaaatgtgtctttcaccgaaactgaaa
ccgaagctcgctgcttctgcgaactacgaacctaacagctgggactatgattacctgctgtcctc
cgacacggacgagtcctcgaagtatacaaaagcaaaagcgaaaaagctggaagccgaagttcgt
cgcgagattaataacgaaaaagcagaatttctgacctgctggaactgattgacaacgtccagc
gcctgggcctgggttaccgttctgagctctgatccgtggtgcgctggatcgcttcgtttcctc
cggcggcttcgatgcggttaaccaagacttcctgcacggtacggcactgtcttccgtctgctg
cgtcaacacgggtttgaggtttctcaggaagcgttcagcggcttcaaagaccaaaccggcaact

Figure 20C


Nucleotide sequence of P. alba pET24a

tcctggagaacctgaaggaagatatcaaagctatcctgagcctgtacgaggccagcttcctggc
tctggaaggcgaaaacatcctggacgagcggaagggttttcgcaatctctcatctgaaagaactg
tctgaagaaaagatcggtaaagagctggcagaacaggtgaaccatgcactggaactgccactgc
atgcgcgtactcagcgtctggaagcagtatggtctatcgaggcctaccgtaaaaaggaggacgc
gaatcaggttctgctggagctggcaattctggattacaacatgatccagtctgtataccagcgt
gatctgcgtgaaacgtcccgttgggtggcgtcgtgtgggtctggcgaccaaactgcactttgctc
gtgaccgcctgattgagagcttctactgggccgtgggtgtagcattcgaaccgcaatactccga
ctgcctgaactccgtcgcaaaaatgttttcttctcgtaaccattatcgacgatatctacgatgta
tacggcaccctggacgaactggagctgtttactgatgcagttgagcgttgggacgtaaacgcca
tcaacgacctgccggattacatgaaactgtgctttctggctctgtataacactattaacgaaat
cgcctacgacaacctgaaagataaagggtgagaacatcctgccgtatctgaccaaagcctgggct
gacctgtgcaacgctttcctgcaagaagccaagtggctgtacaacaaatctactccgaccttg
acgaactacttcggcaacgcatggaaatcctcttctggcccgtgcaactggtgttcgcttactt
cgctgtcgtgcagaacattaaaaaggaagagatcgaaaacctgcaaaaataccatgacaccatc
tctcgtccttcccatatcttccgtctgtgcaatgacctggctagcgcgtctgcggaaattgcgc
gtggtgaaaccgcaaatagcgtttcttgttacatgcgcactaaaggatctccgaagaactggc
taccgaaagcgtgatgaatctgatcgatgaaacctggaaaaagatgaacaaggaaaaactgggt
ggtagcctgttcgcgaaaccgttcgtggaaaccgcgatcaacctggcacgtcaatctcactgca
cttatcataacggcgacgcgcataacctctccggatgagctgacccgcaaacgcgttctgtctgt
aatcactgaaccgattctgccgtttgaacgctaaggatccgaattcgagctccgtcgacaagct
tgcggccgcactcgagcaccaccaccaccactgagatccggtgctaacaaagcccgaaag
gaagctgagttggctgctgccaccgctgagcaataactagcataaacccttggggcctctaaac
gggctcttgaggggttttttgcgtgaaaggaggaactatatccggat

Figure 21

10% NUPAGE Bis-Tris Gel, Reduced

Lane	Gel No. 1 Volume	Description
1	10	MW Marker
2	1.5	Purif IsoS 0.5 ug
3	3.1	Purif IsoS 1 ug
4	6.3	Purif IsoS 2 ug
5	15	P alba Uninduced
6	15	P alba Induced
7	15	P alba Lysate
8	15	P alba Supernatant
9	10	P alba Pellet
10	0	Blank
Function		Coomassie Stain



SUBSTITUTE SHEET (RULE 26)




Figure 23




Figure 24




Figure 25

Amino acid sequence of "MEA" variant of *P. alba* IspS (SEQ ID NO:122)

MEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLELIDNVQ
RLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFVVSQEAFTSGFKDQNGN
FLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVNHAELEPL
HRRRTQRLEAVWSIEAYRKKEDANQVLELAILDYNMIQSVYQRDLRETSRWRRRVGLATKLHFA
RDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLELELFTDAVERWDVNA
INDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNKSTPTF
DDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEEENLQKYHDTISRPSHIFRLCNDLASASAEIA
RGETANSVSCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLARQSHC
TYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 26A

Nucleotide sequence of pDu39 (SEQ ID NO:123)

tggcgaatgggacgcgcctgt agcggcgcat taagcgcggcgggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgcctagcgcgcgcctcttctgctttcttcccttcttctctcg
ccacgttcgcccggctttccccgtcaagctctaaatcgggggctccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaaacttgattagggatgatggttcacgtagtgggcatcg
ccctgatagacgggttttctgccccttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttggt
tatttttctaaatacattcaaataatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaataagg
tatcaagtgaagaaatcaccatgagtgacgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatcaacc
aaaccgttattcattcgtgattgcgcctgagcgcgagcgaataacgcgatcgctgttaaaaggac
aattacaaacaggaatcgaatgcaaccggcgcaggaacactgccagcgcacacaataattttc
acctgaatcaggatattcttctaataacctggaatgctgttttcccggggatcgagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaatccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcacatgggcttcccatacaatcgatagattgtcgacactgattgcccagaca
ttatcgcgagcccatttatacccatataaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttcccggttgaaatggctcataacaccccttgattactgtttatgtaagcaga
cagttttattgttcattgacaaaaatcccttaacgtgagttttcggttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggtaactggcttcagcagagcgcagataccaaataactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgcctacatacctcgctctgctaatacctgttacca
gtggctgctgccagtggcgataagtctgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgcagcggctggggtgaacggggggttcgtgcacacagcccagcttggaagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaaagggaga
aaggcggacaggtatccggttaagcggcaggggtcggaacaggagagcgcacaggggagcttccag
ggggaacgcctggtatctttatagtctgtcggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcaggggggaggagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttctggccttttctcacatgttctttcctgcgttatcccctgatctctgtgg
ataaccgtattaccgcctttgagtgcgtgataccgctcgccgcagccgaacgaccgagcgcag
cgagtgcgtgagcgcaggaagcgcgaagagcgcctgatgcgggtattttctccttacgcatctgtgc
ggtatttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacacccgccaacac
ccgctgacgcgcctgacgggcttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgaattcacagatgtctgcctgttcatccgcgtccagct
cgttgagtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaaggcggt
tttttctgtttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata

Figure 26B

Nucleotide sequence of pDu39

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacaaactggcggatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggtgtccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcgcaggtcgcagacgttttgcagcagcagtcgcttca
cgttcgctcgcgatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgcccggcgataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtgacgaaggcttgagcgagggcggtgcaagat
tccgaataaccgaagcgacaggccgcatcatcgctcgcgctccagcgaaagcggtcctcgccgaaa
atgacccagagcgctgcccgcacctgtcctacgagttgcatgataaagaagacagtcataagt
cgggcagcatagtcatgccccgcgcccaccggaaggagctgactgggttgaaaggctctcaagg
catcggtcgagatcccggtgcctaatgagtgagctaacttacattaattgcgttgcgctcactg
cccgctttccagtcgggaaacctgtcggtccagctgcattaatgaatcgccaacgcgcgggga
gaggcggtttgcgtattgggcgccagggtggtttttctttccaccagtgcagcgggcaacagct
gattgcccttcaccgcctggcctgagagagttgcagcaagcggtccacgctgggttgccccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatccgcaccaacgcgcagcccggactcggtaatggcgcgcattg
cgcccagcgccatctgatcggtggcaaccagcatcgagtggaacgatgccctcattcagcat
ttgcatgggtttgttgaaaaccggacatggcactccagtcgccttccggtccgctatcggtga
atttgattgcgagtgagatatttatgccagccagccagacgcagacgcgcggagacagaactta
atgggcccgcctaacagcgcatgttgcgtggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaataactgttgatgggtgtctgggtcagagacatcaagaaat
aacgcccgaacattagtgaggcagcttccacagcaatggcatcctgggtcatccagcggatagt
taatgatcagccactgacgcgttgccgcgagaagattgtgcaccgcccgtttacaggettgcac
gccgcttcgttctaccatcgacaccaccagctggcaccagttgatcggcgcgagatttaate
gccgcgacaatttgcgacggcgctgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgcgcacgcggttgggaatgtaattcagctccgccatcgccgc
ttccactttttcccgcttttcgcagaaacgtggctggcctgggtccaccacgcgggaaacggtc
tgataagagacaccggcataactctgcgacatcgataaacgttactgggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataaccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggt
gaggccgttgagcaccgcccgcgcaagggaatggtgcatagcaaggagatggcgcccaacagtc
ccggccacggggcctgccaccataccacgcgaaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaatataacgac
tcactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatggaagctcgtcggtctgcgaactacgaacctaacagctgggactatg
attacctgctgtcctccgacacggacgagtcctcgaagtatacaaa'gacaaagcgaaaaagct
ggaagccgaagttcgtcgcgagattaataacgaaaaagcagaatttctgaccctgctggaactg
attgacaacgtccagcgctgggctgggttaccgtttcgagctctgatccgtgggtgcgctgg
atcgcttcgttctccggcggttcgatgcggttaaccaagacttccctgcacggtacggcact
gtcttccgtctgctgcgtcaacacggttttgaggtttctcaggaagcggttcagcggcttcaa
gacaaaacggcaacttccctggagaacctgaagggaagatatcaaagctatcctgagcctgtacg

Figure 26C

Nucleotide sequence of pDu39

aggccagcttcctggctctggaaggcgaaaacatcctggacgaggcgaaaggttttcgcaatctc
tcatctgaaagaactgtctgaagaaaagatcggtaaaagagctggcagaacaggtgaaccatgca
ctggaactgccactgcatcgccgtactcagcgtctggaagcagtatggctctatcgaggcctacc
gtaaaaaggaggacgcgaatcaggttctgctggagctggcaattctggattacaacatgatcca
gtctgtataccagcgtgatctgctgaaacgtcccgctgggtggcgtcgtgtgggtctggcgacc
aaactgcactttgctcgtgaccgcctgattgagagcttctactgggccgtgggtgtagcattcg
aaccgcaatactccgactgccgtaactccgtcgcaaaaatgttttcttcgtaaccattatcga
cgatatctacgatgtatacggcacccctggacgaactggagctgtttactgatgcagttgagcgt
tgggacgtaaacgccatcaacgacctgccggattacatgaaactgtgctttctggctctgtata
acactattaacgaaatcgctacgacaacctgaaagataaagggtgagaacatcctgccgtatct
gaccaaagcctgggctgacctgtgcaacgctttcctgcaagaagccaagtggctgtacaacaaa
tctactccgaccttgacgactacttcggcaacgcacggaaatcctcttctggcccgtgcaac
tgggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatcgaaaacctgcaaaa
ataccatgacaccatctctcgtccttcccatatcttccgtctgtgcaatgacctggctagcgcg
tctgcggaaattgcgcggtggtgaaaccgcaaatagcgtttcttggtacatgcgcactaaaggta
tctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacctggaaaaagatgaa
caaggaaaaactgggtggtagcctgttcgcgaaaccgttcgtggaaaccgcgatcaacctggca
cgtcaatctcactgcacttatcataacggcgacgcgcataacctctccggatgagctgacccgca
aacgcgttctgtctgtaatactgaaccgattctgccgtttgaacgctaaggatccgaattcga
gctccgtcgacaagcttgcgccgcactcgagcaccaccaccaccactgagatccggtgc
taacaaagcccgaaggaagctgagttggctgctgccaccgctgagcaataactagcataaacc
cttggggcctctaaacgggtcttgaggggttttttgctgaaaggaggaactatatccggat

Figure 27

Amino acid sequence of truncated "MSV" variant of *P. alba* IspS (SEQ ID NO:124)

MSVSTENVSFTEETEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEK
AEFLTLLLELIDNVQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFV
SQEAFSGFKDQNGNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGK
ELAEQVNHALELPLHRRRTQRLEAVWSIEAYRKEDANQVLELAILDYNNMIQSVYQDLRETSR
WRRRVGLATKLHFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLDL
ELFTDAVERWDVNAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFL
QEAKWLYNKSTPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEEENLQKYHDTISRPSHIF
RLCNDLASASAEIARGETANSVSCYMRKGISEELATESVMNLIDETWKKMNKEKLGGSFPAKP
FVETAINLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 28A

Nucleotide sequence of pDu41 (pET24a-P.alba (MSV) Untagged) (SEQ ID NO:125)

tggcgaatgggacgcgcctgtagcggcgcattaagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccctagcggcgctcctttcgctttcttcccttctttctcg
ccacgttcgccggctttccccgtcaagctctaaatcgggggctccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaaacttgattaggtgatggttcacgtagtgggcatcg
ccctgatagacggtttttcgccctttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaatatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaataagg
tatcaagtgagaaatcaccatgagtgcgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatcaacc
aaaccgttattcattcgtgattgcgcctgagcgagacgaaatacgcgatcgctgttaaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatcaacaatattttc
acctgaatcaggatattcttctaataacctggaatgctgttttcccggggatcgagtggtgagt
aaccatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatcgggcttcccatacaatcgatagattgtcgacactgattgcccagaca
ttatcgcgagccatttataccatataaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttcccggtgaatatggctcataacaccccttgtattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcggttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgcagataccaaatactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgctacatacctcgctctgctaatacctgttacca
gtggctgctgccagtggcgataaagtcgtgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgcgagcgtcgggctgaacggggggttcgtgcacacagcccagcttggaagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaaagcgccacgcttcccgaagggaga
aaggcgagacaggtatccggtgaacggcgagggtcggaacaggagagcgacaggggagcttccag
ggggaacgcctggtatctttatagtctgtcggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatcccctgatctctgtgg
ataaccgtattaccgcctttgagtgcgctgataccgctcgccgcagccgaacgaccgagcgag
cgagtgcgtgagcgaggaagcggaagagcgctgatgcggtattttctccttacgcatctgtgc
ggtatttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagttaa
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacacccgccaacac
ccgctgacgcgcctgacgggcttgcgtgctcccgcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcattcacagatgtctgcctgttcatccgcgtccagct
cggtgagtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaagggcggt

Figure 28B

Nucleotide sequence of pDu41 (pET24a-P.alba (MSV) Untagged)

tttttcctggttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata
ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccgggttactgg
aacgttgtgagggtaaacactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtaggtgttccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcaggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgtgtgctcaggtcgcagacgttttgcagcagcagtcgcttca
cggttcgctcgcgtatcgggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgatcatgcgcacccgtggggccgcatgcccggcgataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtgacgaaggcttgagcagggcggtgcaagat
tccgaataaccgcaagcgacaggccgatcatcgtcgcgctccagcgaaagcggtcctcgccgaaa
atgacccagagcgctgcgggcacctgtcctacgagttgcatgataaagaagacagtcataagt
cgggcagcatagtcatgccccgcgcccaccggaaggagctgactgggttgaaggctctcaagg
catcggtcgagatcccggtgcctaatgagtgcgtaacttacattaattgcgttgcgctcactg
cccgtttccagtcgggaaacctgtcgtgccagctgcattaatgaatcgccaacgcgcgggga
gaggcggtttgcgtattggcgccagggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttcaccgcctggcctgagagagttgcagcaagcggtccacgctggtttgccccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcgggtatcg
tcgtatcccactaccgagatatccgcaccaacgcgcagcccggactcggtaatggcgcgcat
cgcccagcgccatctgatcgttggcaaccagcatcgcagtggaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttccggttccgctatcggtga
atltgatgagagtgagatatltatgccagccagccagacgcagacgcgcgagacagaactta
atgggcccgttaacagcgcgatttgcgtggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaataactgttgatgggtgtctggtcagagacatcaagaaat
aacgcccgaacattagtgaggcagcttccacagcaatggcatcctggtcatccagcggatagt
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggcttcgac
gccgttctgttctaccatcgacaccaccagctggcaccagttgatcggcgcgagatttaac
gccgcgacaatttgcgacggcgcggtgcaggggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgtgtgccacgcggttgggaatgtaattcagctccgccatcgccgc
ttccactttttccgcggttttcgcagaaacgtggctggcctggttcaccacgcgggaaacgggtc
tgataagagacaccggcataactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataaccgcgaaagggttttgccgcatcagatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcccgcgcaagggaatggtgcattgcaaggagatggcgcccaacagctccc
ccggccacggggcctgccaccatacccacgcggaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatcagac
tactataggggaattgtgagcggataacaattccccctctagaaataattttgtttaactttaa
gaaggagatatacatatgagcgtgtccaccgaaaatgtgtctttcaccgaaactgaaaccgaag
ctcgtcgttctgcgaactacgaacctaaccagctgggactatgattacctgctgtcctccgacac
ggacgagtcctatcgaagtatacaaaagacaaagcgaaaaagctggaagccgaagtctgtcgcgag
attaataacgaaaaagcagaatttctgacctgctggaactgattgacaacgtccagcgcctgg
gcctgggttaccgttctgagtcctgatatccgtggtgcgctggatcgcttcgtttcctccggcgg
cttcgatgcggtaaccaagacttcctgcacggtacggcactgtcttccgctctgctgcgtcaa

Figure 28C

Nucleotide sequence of pDu41 (pET24a-P.alba (MSV) Untagged)

cacgggttttgagggtttctcaggaagcggttcagcggcttcaaagacccaaacggcaacttcctgg
agaacctgaaggaagatatcaaagctatcctgagcctgtacgaggccagcttcctggctctgga
aggcgaaaacatcctggacgaggcgaagggttttcgcaatctctcatctgaaagaactgtctgaa
gaaaagatcggtaaagagctggcagaacaggtgaaccatgcactggaactgccactgcatcgcc
gtactcagcgtctggaagcagtatggtctatcgaggcctaccgtaaaaaggaggacgcgaatca
ggttctgctggagctggcaattctggattacaacatgatccagtctgtataccagcgtgatctg
cgtgaaacgtcccgttgggtggcgtcgtgtgggtctggcgaccaaactgcactttgctcgtgacc
gcctgattgagagcttctactgggccgtgggtgtagcattcgaaccgcaataactccgactgccg
taactccgtcgcaaaaatgttttcttcgttaaccattatcgacgatatctacgatgtatacggc
accctggacgaactggagctgttactgatgcagttgagcgttgggacgtaaacgccatcaacg
acctgccggattacatgaaactgtgctttctggctctgtataacactattaacgaaatcgcta
cgacaacctgaaagataaagggtgagaacatcctgccgtatctgaccaaagcctgggctgacctg
tgcaacgctttcctgcaagaagccaagtggctgtacaacaaatctactccgacctttgacgact
acttcggcaacgcattgaaatcctcttctggcccgtgcaactgggtgttcgcttacttcgctgt
cgtgcagaacattaaaaaggaagagatcgaaaacctgcaaaaataccatgacaccatctctcgt
ccttcccatatcttccgtctgtgcaatgacctggctagcgcgtctgcggaaattgcgcgtgggtg
aaaccgcaaatagcgtttcttgttacatgcgcactaaagggtatctccgaagaactggctaccga
aagcgtgatgaatctgatcgatgaaacctggaaaaagatgaacaaggaaaaactgggtggtagc
ctgttcgcgaaaccgttcgtggaaaccgcgatcaacctggcacgtcaatctcactgcacttatc
ataacggcgacgcgcatacctctccggatgagctgacccgcaaacgcgttctgtctgtaatcac
tgaaccgattctgccgtttgaacgctaaggatccgaattcgagctccgtcgacaagcttgccgc
cgcactcgagcaccaccaccaccactgagatccggctgctaacaagcccgaagggaagct
gagttggctgctgccaccgctgagcaataactagcataaacccttggggcctctaaacgggtct
tgaggggttttttgctgaaaggaggaactatatccggat

Figure 29

Amino acid sequence of truncated "MVS" Variant of *P. alba* IspS (SEQ ID NO:126)

MVSTENVSFTEETETEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKA
EFLTLLLELIDNVQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFVS
QEA FSGFKDQNGNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKE
LAEQVNHALELPLHRRRTQRLEAVWSIEAYRKEDANQVLELAILDYNMIQSVYQRDLRETSRW
WRRVGLATKLHFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTI IDDIYDVYGTLDLE
LFTDAVERWDVNAINDLDPYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQ
EAKWLYNKSTPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEE IENLQKYHDTISRPSHIFR
LCNDLASASAEIARGETANSVSCYMRKGI SEELATESVMNLI DETWKKMNKEKLGGSLFAKPF
VETAINLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 30A

Nucleotide sequence of pDu43 (pET24a-P.alba (MVS) Untagged) (SEQ ID NO:127)

tggcgaatgggacgcgcctgtagcggcgcat taagcgcggcggggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccttagcgcgcgtcctttcgctttcttcccttcccttctcg
ccacgttcgcgcgtttccccgtcaagctctaaatcgggggctcccttaggggtccgatttag
tgctttacggcacctcgaccccaaaaaacttgattaggggtgatgggtcacgtagtgggccatcg
ccctgatagacgggttttcgccccttgacgttggagtcacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggattttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaatatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttccctcgtcaaaaaataaggt
tatcaagtgagaaatcaccatgagtgacgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcataacc
aaaccgttattcattcgtgattgcgcctgagcgcagacgaaatacgcgatcgtgtttaaaggac
aattacaaacaggaatcgaatgcaaccggcgcaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaatacctggaatgctgttttccgggggatcgcagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatacgggttcccatacaatcgatagattgtgcacactgattgcccagaca
ttatcgcgagccatttataccatataaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttccggttgaatatggctcataaacaccccttgtattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcggttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgttggtttgtttgccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgcagataccaaataactgtccttctagtgtagecgtagtt
aggccaccacttcaagaactctgtagcaccgcctacatacctcgtctctgctaatectgttacca
gtggctgctgccagtggcgataagtcgtgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgcagcggtcgggtgaacgggggggttcgtgcacacagcccagcttggagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgcacagcttcccgaagggaga
aaggcggacaggtatccggtgaacggcaggggtcggaacaggagagcgcacagaggagcttccag
ggggaacgcctggtatctttatagtcctgtcgggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcagggggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatccctgattctgtgg
ataaccgtattaccgcctttgagtgcgctgataccgctcgcgcgagccgaacgaccgagcgcag
cgagtgcgtgagcgcggaagcgaagagcgcctgatgcggtattttctccttacgcatctgtgc
ggtatttcacaccgcatataggtgcactctcagtaaatctgctctgatgccgcatagttaag
ccagtataactccgctatcgctacgtgactgggtcatggctgcgccccgacaccgcgaacac
ccgctgacgcgcctgacgggttgtctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cgttgagtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaaggggcgt
ttttcctgtttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata

Figure 30B

Nucleotide sequence of pDu43 (pET24a-P.alba (MVS) Untagged)

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacaactggcggatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggtgtccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttggtgctcaggtcgcagacgttttgcagcagcagtcgcttca
cgttcgctcgcgatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgcccggcgaataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtgacgaaggcttgagcgaaggcggtgcaagat
tccgaataccgcaagcgacaggccgcatcatcgtcgcgctccagcgaaagcggtcctcgccgaaa
atgacccagagcgctgcccggcacctgtcctacgagttgcatgataaagaagacagtcataagt
cggcgacgatagtcgccccgcgcccaccggaaggagctgactgggttgaaaggctctcaagg
catcggtcgagatcccgggtgcctaatagtgagctaaactacattaattgcgttgcgctcactg
cccgctttccagtcgggaaacctgtcgtgccagctgcattaatgaatcggccaaacgcgcgggga
gaggcggtttgcgtattgggcgccagggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttcaccgcctggccctgagagagttgcagcaagcggtccacgctgggttgcccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatatccgcaccaacgcgcagcccggactcggtaatggcgcgcattg
cgcccagcgccatctgatcgttggaaccagcatcgagtgggaaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttcgctatcggtga
atttgattgcgagtgagatatattatgccagccagccagacgcagacgcgcgcgagacagaactta
atgggcccgcctaacagcgcgatttgctggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaaataactgttgatgggtgtctgggtcagagacatcaagaaat
aacgccggaacattagtgcaggcagcttccacagcaatggcatcctgggtcatccagcggatagt
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggttcgcac
gccgcttcggttctaccatcgacaccaccacgctggcaccagttgatcggcgcgagatttaate
gccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgcccaatcagcaacg
actgtttgccgcaggttggtgtgccacgcggttggaatgtaattcagctccgccatcgccgc
ttccactttttccgcgttttcgcagaaacgtggctggcctgggtcaccacgcgggaaacgggtc
tgataagagacaccggcataactctgcgacatcgtataacgttactgggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcccggcaaggaatggtgcatgcaaggagatggcgcccaacagtcctc
ccggccacggggcctgccaccatacccacgcgcgaaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatacgac
tcactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatggtgtccaccgaaaaatgtgtctttcaccgaaactgaaaccgaagctc
gtcgttctgcgaactacgaacctaacagctgggactatgattacctgctgtcctccgacacgga
cgagtcctcgaagtatacaaagacaaagcgaaaaagctggaagccgaagttcgtcgcgagatt
aataacgaaaaagcagaatttctgaccctgctggaactgattgacaacgtccagcgctgggccc
tgggttacggttcgagctcgatatccgtgggtgcgctggatcgcttcgtttcctccggcggtt
cgatgcggttaaccaagacttccctgcacggtagggcactgtcttccgtctgctgcgtcaacac
ggttttgaggtttctcaggaagcggttcagcggttcaaagaccaaaccggcaacttccctggaga

Figure 30C

Nucleotide sequence of pDu43 (pET24a-P.alba (MVS) Untagged)

acctgaaggaagatatcaaagctatcctgagcctgtacgaggccagcttcctggctctggaagg
cgaaaacatcctggacgaggcgaaggttttcgcaatctctcatctgaaagaactgtctgaagaa
aagatcggtaaaagagctggcagaacaggtgaaccatgcactggaactgccactgcatcgccgta
ctcagcgtctggaagcagtatggctctatcgaggcctaccgtaaaaaggaggacgcgaatcaggt
tctgctggagctggcaattctggattacaacatgatccagtctgtataccagcgtgatctgcgt
gaaacgtcccgttgggtggcgtcgtgtgggtctggcgaccaaactgcactttgctcgtgaccgcc
tgattgagagcttctactgggccgtgggtgtagcattcgaaccgcaataactccgactgccgtaa
ctccgtcgcaaaaatgttttcttctcgtaaccattatcgacgatatctacgatgtatacggcacc
ctggacgaactggagctgttactgatgcagttgagcgttgggacgtaaacgccatcaacgacc
tgccggattacatgaaactgtgcttcttggtctgtataaactattaacgaaatcgctacga
caacctgaaagataaaggtgagaacatcctgccgtatctgaccaaagcctgggctgacctgtgc
aacgctttcctgcaagaagccaagtggctgtacaacaaatctactccgaccttgacgactact
tcggcaacgcatggaaatcctcttctggcccgtgcaactgggtgttcgcttacttcgctgtcgt
gcagaacattaaaaaggaagagatcgaaaacctgcaaaaataccatgacaccatctctcgtcct
tcccatatcttccgtctgtgcaatgacctggctagcgcgtctgcggaaattgcgcgtggtgaaa
ccgcaaatagcgtttcttgttacatgcgcactaaagggtatctccgaagaactggctaccgaaag
cgtgatgaatctgatcgatgaaacctggaaaaagatgaacaaggaaaaactgggtggtagcctg
ttcgcgaaaccgttcgtggaaaccgcgatcaacctggcacgtcaatctcactgcacttatcata
acggcgacgcgcataacctctccggatgagctgacccgcaaacgcgttctgtctgtaatactga
accgattctgccgtttgaaacgctaaggatccgaattcgagctccgtcgacaagcttgccggccgc
actcgagcaccaccaccaccactgagatccggctgctaacaaagcccgaaaggaagctgag
ttggctgctgccaccgctgagcaataactagcataaacccttggggcctctaaacgggtcttga
ggggttttttgcgtgaaaggaggaactatatccggat

Figure 31

Amino acid sequence of truncated "MTE" variant of *P. alba* IspS (SEQ ID NO:128)

MTENVSFTEETETARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEF
LTLLLELIDNVQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFEVSQE
AFSGFKDQNGNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELA
EQVNHALELPLHRRRTQRLEAVWSIEAYRKKEDANQVLELAILDYNMIQSVYQDLRETSRWWR
RVGLATKLHFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTI IDDIYDVYGTLDLELELF
TDAVERWDVNAINDLDPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEA
KWLYNKSTPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEEIEENLQKYHDTISRP SHIFRLC
NDLASASAEIARGETANSVSCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVE
TAINLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 32A

Nucleotide sequence of pDu42 (pET24a-P.alba (MTE) Untagged) (SEQ ID NO:129)

tggcgaatgggacgcgcctgtagcggcgcatataagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccctagcggccgctcctttcgctttcttcccttcccttctcg
ccacgttcgcccgtttccccgtcaagctctaaatcgggggtccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaacttgattagggtgatggttcacgtagtgggcatcg
ccctgatagacggtttttcgccctttgacgttgaggtccacgttctttaatagtggactctgt
tccaaactggaacaactcaacctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttggt
tatttttctaaatacattcaaatatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
ttctgtaatgaaggagaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttccctcgtcaaaaaataagg
tatcaagtgaagaatcaccatgagtgcgactgaatccggtgagaatggcaaaagttaatgc
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcataacc
aaaccgttattcattcgtgattgcgcctgagcgcgagcaataacgcgatcgtgtttaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaataacctggaatgctgttttccgggggatcgcagtggtgag
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatcgggcttcccatacaatcgatagattgtcgcacctgattgcccgaca
ttatcgcgagccatttataccatataaatcagcatccatgttggaatttaatecgggcctag
agcaagacgtttccggttgaaatgggtcataaacaccccttgattactgtttatgtaagcaga
cagttttattgttcattgacccaaatccctaacgtgagttttcgttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggtaactggcttcagcagagcgcagataccaaatactgtccttctagtgtagecgtagtt
aggccaccacttcaagaactctgtagcaccgcctacatacctcgtctctgctaatectgttacca
gtggctgctgccagtggcgataagtctgtcttaccgggttggtactcaagacgatagttaccgg
ataaggcgcagcggcgggtgaacggggggttcgtgcacacagcccagcttgagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttccgaagggaga
aaggcggacaggtatccggtgaacggcgagggtcggaacaggagagcgcacgagggagcttccag
ggggaacgcctggtatctttatagtcctgtcgggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatccctgattctgtgg
ataaccgtattaccgcctttgagtgcgctgataccgctcgcgcgagccgaacgaccgagcgcag
cgagtgcgtgagcgcggaagcgaagagcgccctgatgcggatattttctccttacgcatactgtgc
ggtatttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagttaa
ccagatacactccgctatcgctacgtgactgggtcatgggtgcgcccgcacaccgcgaacac
ccgctgacgcgcctgacgggcttgctgctccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagagggttttcacgcgtcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgaattcacagatgtctgcctgttcatccgcgtccagct
cgttgagtttctccagaagcgttaaatgtctggctctcgataaagcggggccatgttaaggcggt
tttttctgtttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata

Figure 32B

Nucleotide sequence of pDu42 (pET24a-P.alba (MTE) Untagged)

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacactggcggatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggtgtccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcaggtcgcagacgttttgcagcagcagtcgcttca
cggttcgctcgcgtatcgggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgccggcgataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtacgaaggcttgagcgaggggcgtaagat
tccgaataccgcaagcgacaggccgcatcatcgctcgcgctccagcgaaagcggtcctcgccgaaa
atgacccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagt
cggcgacgatagtcatgccccgcgcccaccggaaggagctgactgggttgaaggctctcaagg
catcggtcgagatcccggtgcctaatagtgagctaaactacattaattgcgttgcgctcactg
ccgcgtttccagtcgggaaacctgtcgtgccagctgcattaatgaatcggccaaacgcgcgggga
gaggcggtttgcgtattgggcgccagggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttcaccgcctggccctgagagagttgcagcaagcggtccacgctggtttgccccag
caggcgaaaatcctggttgatggtggttaacggcgggataatacatgagctgtcttcggtatcg
tcgtatcccactaccgagatatccgcaccaacgcgcagcccggactcggtaatggcgcgcattg
cgcccagcgccatctgatcgttggcaaccagcatcgagtggaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atattgattgcgagtgagatatttatgccagccagccagacgcagacgcgcccagacagaaacta
atgggcccgcataacagcgcgatttgcgtggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaataatactgttgatgggtgtctggtcagagacatcaagaaat
aacgcgggaacattagtgaggcagcttccacagcaatggcatcctggtcatccagcggatagtc
taatgatcagccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggcttcgac
gccgcttcggttctaccatcgacaccaccacgctggcaccagttgatcggcgcgagatttaatac
gccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgcccaatcagcaacg
actgtttgcccgccagttgttgcgcacgcggttgggaatgtaattcagctccgccatcgccgc
ttccactttttcccgcggttttcgcagaaacgtggctggcctggttaccacgcgggaaacgggtc
tgataagagacaccggcataactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttcggggcgctatcatgccataccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcgcgcgcaagggaatggtgcatgcaaggagatggcgcccaacagtcct
ccggccacggggcctgccaccataccacgcgcaaaacagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatacgac
tactataggggaattgtgagcggataacaattccctctagaaataattttgtttaactttaa
gaaggagatatacatatgaccgaaaatgtgtctttaccgaaactgaaaccgaagctcgtcggtt
ctgcgaactacgaacctaacagctgggactatgattacctgctgtcctccgacacggacgagtc
catcgaagtatacaagacaaagcgaaaaagctggaagccgaagttcgtcgcgagatttaataac
gaaaaagcagaatttctgacctgctggaactgattgacaacgtccagcgccctgggcctgggtt
accgtttcgagctctgatatccgtggtgcgctggatcgcttcgtttcctccggcggttcgatgc
ggttaaccaagacttccctgcacggtacggcactgtcttccgtctgctgcgtcaacacgggtttt
gaggtttctcaggaagcgttcagcggcttcaaagaccaaaccggcaacttctggagaacctga

Figure 32C

Nucleotide sequence of pDu42 (pET24a-P.alba (MTE) Untagged)

aggaagatatcaaagctatcctgagcctgtacgaggccagcttcttggctctggaaggcgaaaa
catcctggacgaggcggaaggttttcgcaatctctcatctgaaagaactgtctgaagaaaagatc
ggtaaagagctggcagaacaggtgaaccatgcactggaactgccactgcatcgccgtactcagc
gtctggaagcagtatggtctatcgaggcctaccgtaaaaaggaggacgcgaatcaggttctgct
ggagctggcaattctggattacaacatgatccagtctgtataccagcgtgatctgctgaaacg
tcccgttgggtggcgtcgtgtgggtctggcgaccaaactgcactttgctcgtgaccgcctgattg
agagcttctactgggccgtgggtgtagcattcgaaccgcaataactccgactgccgtaactccgt
cgcaaaaatgttttctttcgttaaccattatcgacgatatctacgatgtatacggcacccctggac
gaactggagctgtttactgatgcagttgagcgttgggacgtaaacgccatcaacgacctgccgg
attacatgaaactgtgctttctggctctgtataacactattaacgaaatcgctacgacaacct
gaaagataaaggtgagaacatcctgccgtatctgaccaaagcctgggctgacctgtgcaacgct
ttcctgcaagaagccaagtggctgtacaacaaatctactccgacctttgacgactacttcggca
acgcatggaaatcctcttctggcccgtgcaactgggtgttcgcttacttcgctgtcgtgcagaa
cattaaaaagggaagagatcgaaaacctgcaaaaataccatgacaccatctctcgtccttcccat
atcttccgtctgtgcaatgacctggctagcgcgtctgcggaaattgcgcgtgggtgaaaccgcaa
atagcgtttcttgttacatgcgcactaaaggtatctccgaagaactggctaccgaaagcgtgat
gaatctgatcgatgaaacctggaaaaagatgaacaaggaaaaactgggtggtagcctgttcgcg
aaaccgttcgtggaaaccgcgatcaacctggcacgtcaatctcactgcacttatcataacggcg
acgcgcatacctctccggatgagctgacccgcaaaccgcgttctgtctgtaatactgaaccgat
tctgccgtttgaacgctaaggatccgaattcgagctccgtcgacaagcttgcggccgcactcga
gcaccaccaccaccaccactgagatccggctgctaacaaagcccgaagggaagctgagttggct
gctgccaccgctgagcaataactagcataaaccccttggggcctctaaacgggtcttgaggggtt
ttttgctgaaaggagggaactatatccggat

Figure 33

Amino acid Sequence of truncated "MNV" variant of *P. alba* IspS (SEQ ID NO:130)

MNVSFTETETEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLT
LLELIDNVQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFEVSQEAF
SGFKDQNGNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQ
VNHALELPLHRRTRQLEAVWSIEAYRKEDANQVLELAILDYNNMIQSVYQDRLRETSRWRRRV
GLATKLHFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLDLELFTD
AVERWDVNAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKW
LYNKSTPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEEIEENLQKYHDTISRP SHIFRLCND
LASASAEIARGETANSVSCYMRKGI SEELATESVMNLI DETWKKMNKEKLGGSLFAKPFVETA
INLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 34A

Nucleotide sequence of pDu40 (pET24a-P.alba (MNV) Untagged) (SEQ ID NO:131)

tggcgaatgggacgcgcctgtagcggcgcattaagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccttagcgcgcctcctttcgctttcttcccttctttctcg
ccacgttcgcggctttccccgtcaagctctaaatcgggggtccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaaacttgattagggtgatggttcacgtagtgggccatcg
ccctgatagacggtttttcgccctttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggattttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaataatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaataaggt
tatcaagtgagaaatcaccatgagtgcgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcataacc
aaaccgttattcattcgtgattgcgcctgagcgcgagacgaaatacgcgatcgtgtttaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaataacctggaatgctgtttcccggggatcgagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatacgggttcccatacaatcgatagattgtcgcacctgattgcccagaca
ttatcgcgagccatttataccatataaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttcccggttgaaatggctcataaacaccccttgattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcggttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgcagataccaaatactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgcctacataacctcgctctgctaactcctgttacca
gtggctgctgccagtggcgataagtcgtgtcttaccgggttggtactcaagacgatagttaccgg
ataaggcgcagcgtcgggtgaacggggggttcgtgcacacagcccagcttgagcgaacgcac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaagggaga
aaggcggacaggtatccggttaagcgcgaggggtcggaacaggagagcgcacagaggagcttccag
ggggaacgcctggtatctttatagtcctgtcggtttcgccacctctgacttgagcgtcgatt
tttgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatcccctgattctgtgg
ataaccgtattaccgcctttgagtgcgctgataccgctcgccgcagccgaacgaccgagcgcag
cgagtcagtgcgaggaagcgaagagcgcctgatgcggtattttctccttacgcatactgtgc
ggtatttcacaccgcataatggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagatacactccgctatcgctacgtgactgggtcatggctgcgccccgacacccgccaacac
ccgctgacgcgcctgacgggttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgaattcacagatgtctgcctgttcatccgcgtccagct
cggttaggtttctccagaagcgttaatgtctggcttctgataaagcgggccaatgttaaggcggt
tttttctgtttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata

Figure 34B

Nucleotide sequence of pDu40 (pET24a-P.alba (MNV) Untagged)

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccgggttactgg
aacgttgtgagggtaaaacaactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtaggtgtccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcaggtgcgagacgttttgcagcagcagtcgcttca
cgcttcgctcgcgtatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgcccggcgataatggcct
gcttctcgcgaaacgtttggtggcgggaccagtgacgaagggttgagcagagggcggtgcaagat
tccgaataaccgcaagcgacaggccgcatcatcgtcgcgctccagcgaaagcggtcctcgccgaaa
atgaccagagcgctgcccgcacctgtcctacgagttgcatgataaagaagacagtcataagt
cggcgacgatagtcagccccgcgcccaccggaaggagctgactgggttgagggtctcaagg
catcggtcgagatcccggtgctaatgagtgagctaaactacattaattgcgttgcgctcactg
cccgctttccagtcgggaaacctgtcgtgccagctgcattaatgaatcgcccaacgcgcgggga
gaggcggtttgcgtattgggcgccagggtggtttttcttttccacagtgagacgggcaacagct
gattgcccttcaccgctggcctgagagagttgcagcaagcggtccacgctgggttgcccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatataccgcaccaacgcgcagcccggactcggtaatggcgcgcatg
cgcccagcgccatctgatcgttggcaaccagcatcgagtgggaaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atgtgatgagtgagatatttatgccagccagccagacgcagacgcgcgcgagacagaactta
atgggcccgcctaacagcgcgatttgcgtggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgctctcatgggagaaaaataatactgttgatgggtgtctggtcagagacatcaagaaat
aacgcgggaacattagtgagggcagcttccacagcaatggcatcctgggtcatccagcggatagt
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggcttcgac
gccgcttcgcttaccatcgacaccaccacgctggcaccacagttgatcggcgcgagatttaac
gccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgcgcacgcggttgggaatgtaattcagctccgccatcgccgc
ttccactttttcccgcgttttcgcagaaacgtggctggcctggttcaccacgcgggaaacgggtc
tgataagagacaccggcatactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataccgcgaaagggttttgccgcatcagatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggt
gaggccgttgagcaccgcccgcgcaaggaatggtgcatgcaaggagatggcgcccaacagctccc
ccggccaacggggcctgccaccatacccacgcggaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatcagac
tcactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatgaatgtgtctttcaccgaaactgaaaccgaagctcgtcgttctgcga
actacgaacctaacagctgggactatgattacctgctgtcctccgacacggacgagtcctcga
agtatacaaaagacaaagcgaaaaagctggaagccgaagttcgtcgcgagatttaataacgaaaaa
gcagaatttctgaccctgctggaactgattgacaacgtccagcgctgggcctgggttaccgtt
tcgagctcgatatacgtggtgcgctggatcgcttcgtttcctccggcggttcgatgcggtaac
caagacttccctgcacggtacggcactgtctttccgtctgctgcgtcaacacgggttttgaggtt
tctcaggaagcgttcagcggcttcaaagacaaaacgggaacttccctggagaacctgaaggaag

Figure 34C

Nucleotide sequence of pDu40 (pET24a-P.alba (MNV) Untagged)

atatcaaagctatcctgagcctgtacgaggccagcttcctggctctggaaggcgaaaacatcct
ggacgaggcggaaggttttcgcaatctctcatctgaaagaactgtctgaagaaaagatcggtaaa
gagctggcagaacaggtgaaccatgcactggaactgccactgcatcgccgtactcagcgtctgg
aagcagtatggtctatcgaggcctaccgtaaaaaggaggacgcgaatcaggttctgctggagct
ggcaattctggattacaacatgatccagtctgtataccagcgtgatctgctgaaacgtcccgt
tggctggcgtcgtgtgggtctggcgaccaaactgcactttgctcgtgaccgcctgattgagagct
tctactgggcccgtgggtgtagcattcgaaaccgcaataactccgactgccgtaactccgtcgcaaa
aatgttttcttctcgtaaccattatcgacgatatctacgatgtatacggcaccctggacgaactg
gagctgtttactgatgcagttgagcgttgggacgtaaacgccatcaacgacctgccggattaca
tgaactgtgctttctggctctgtataaactattaacgaaatcgctacgacaacctgaaaga
taaaggtgagaacatcctgccgtatctgaccaaagcctgggctgacctgtgcaacgctttcctg
caagaagccaagtggctgtacaacaaatctactccgacctttgacgactacttcggcaacgcat
ggaaatcctcttctggcccgtgcaactgggtgttcgcttacttcgctgtcgtgcagaacattaa
aaaggaagagatcgaaaacctgcaaaaataccatgacaccatctctcgtccttcccatatcttc
cgtctgtgcaatgacctggctagcgcgtctgcggaaattgcgcgtgggtgaaaccgcaaatagcg
tttcttggttacatgcgcactaaaggtatctccgaagaactggctaccgaaagcgtgatgaatct
gatcgatgaaacctggaaaaagatgaacaaggaaaaactgggtggtagcctgttcgcgaaaccg
ttcgtggaaaccgcgatcaacctggcacgtcaatctcactgcacttatcataacggcgacgcgc
atacctctccgatgagctgacccgcaaacgcgttctgtctgtaatactgaaccgatctctgcc
gtttgaacgctaaggatccgaattcgagctccgtcgacaagcttgcgcccgcaactcgagcacca
ccaccaccaccactgagatccggctgctaacaaagcccgaaggaagctgagttggctgctgcc
accgctgagcaataactagcataaacccttggggcctctaaacgggtcttgaggggttttttgc
tgaaggagggaactatatccgat

Figure 35




Figure 36

Amino acid sequence of *P. alba* MEA(+)-TEV (SEQ ID NO:132)

MEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLLELIDNVQ
RLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFVVSQEA FSGFKDQNGN
FLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVNHAELEPL
HRRTQRLEAVWSIEAYRKKEDANQVLELAILDYNMIQSVYQDLRETSRWRRVGLATKLHFA
RDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLELELFTDAVERWDVNA
INDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNKSTPTF
DDYFGNAWKSSSGPLQLVFAYFAVQNIKKEE IENLQKYHDTISRPSHIFRLCNDLASASAEIA
RGETANSVSCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLARQSHC
TYHNGDAHTSPDELTRKRVLSVITEPILPFERENLYFQGLEHHHHHH

Figure 37A

Nucleotide sequence of MD09-163 (pET24a-P. alba MEA(+)-TEV) (SEQ ID NO:133)

tggcgaatgggacgcgcctgtagcggcgcat taagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccctagcgcgcgctcctttcgctttcttcccttccctttctcg
ccacggttcgcccgtttccccgtcaagctctaaatcgggggctccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaacttgattagggtgatggttcacgtagtgggccatcg
ccctgatagacggtttttcgccctttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggattttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacggtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaataatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaaactcaccgaggcaggtccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaaaggt
tatcaagtgagaaatcaccatgagtgcgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcataacc
aaaccgttattcattcgtgattgcgctgagcgcgagcaaatacgcgatcgtgtttaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaataacctggaatgctgtttcccggggatcgcagtggtgagt
aaccatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatacgggttcccatacaatcgatagattgtcgacactgattgcccagaca
ttatcgcgagccattttatacccatataaatcagcatccatgttggaaatttaatcgcggcctag
agcaagacggtttcccggtgaatatggctcataaacacccttgattactgtttatgtaagcaga
cagttttattgttcatgacccaaatcccttaacgtgagttttcggtccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgcagataccaaataactgtccttctagtgtagccgtagt
aggccaccacttcaagaactctgtagcaccgcctacatacctcgtctctgctaactcctgttacca
gtggctgctgccagtgccgataaagtcgtgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgcagcggtcgggctgaacggggggttcgtgcacacagcccagcttgagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgcacagcttcccgaagggaga
aaggcggacaggtatccggttaagcggcagggtcggaacaggagagcgcacgagggagcttccag
ggggaacgcctggtatctttatagtcctgtcggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatcccctgattctgtgg
ataaccgtattaccgcctttgagtgcgctgataccgctcgcgcgcagccgaacgacgagcgcag
cgagtgcgtgagcgaaggaagcgaagagcgcctgatgcgggtattttctccttacgcatctgtgc
ggtattttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagtttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacacccgccaacac
ccgctgacgcgcctgacgggcttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cgttgagtttctccagaagcgttaatgtctggcttctgataaagcgggcatgttaaggcggt

Figure 37B

Nucleotide sequence of MD09-163 (pET24a-P. alba MEA(+)-TEV)

tttttctgtttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata
ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggtgtccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcaggcgctgacttccgcgttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcaggtcgcagacggtttgcagcagcagtcgcttca
cgttcgctcgcgatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgatcatgcgcacccggtggggccgcatgccggcgataatggcct
gcttctcgccgaaacgttgggtggcgggaccagtgacgaaggcttgagcgagggcggtgcaagat
tccgaataaccgcaagcgacaggccgatcatcgctcgcgctccagcgaaagcggtcctcgccgaaa
atgacccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagtg
cggcgacgatagtcatgccccgcgcccaccggaaggagctgactgggttgaaggctctcaaggg
catcggtcgagatcccggtgctaatgagtgagctaacttacattaattgcgttgcgctcactg
cccgtttccagtcgggaaacctgtcggtccagctgcattaatgaatcgccaaacgcgcgggga
gaggcggtttgcgtattggcgccagggtggtttttctttccaccagtgcagcgggcaacagct
gattgcccttcaccgctggcctgagagagttgcagcaagcggtccacgctggtttgccccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatccgcaccaacgcgcagcccggactcggtaatggcgcgattg
cgcccagcgccatctgatcggtggcaaccagcatcgagtggaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atttgattgcgagtgagatatttatgccagccagccagacgcagacgcgcgagacagaactta
atggggccgctaacagcgcgatttgcgtggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaataatactgttgatgggtgtctggtcagagacatcaagaaat
aacgcccgaacattagtgaggcagcttccacagcaatggcatcctggtcatccagcggatagt
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggttcgcac
gccgcttcgttctaccatcgacaccaccacgctggcaccagttgatcggcgcgagatttaate
gccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgcccaatcagcaacg
actgtttgcccgccagttgttgcacgcgggttgggaatgtaattcagctccgccatcgccgc
ttccactttttcccgcgttttcgcagaaacgtggctggcctggttcaccacgcgggaaacgggtc
tgataagagacaccggcataactctgcgacatcgatataacgttactgggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcccgcgcaaggaatggtgcatgcaaggagatggcgcccaacagtcctc
ccggccacggggcctgccaccataccacgcgcgaaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcgcgcatataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaatataacgac
tactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatggaagctcgtcgttctgcgaactacgaacctaacagctgggactatg
attacctgctgtcctccgacacggacgagtcctcgaagtatacaagacaaagcgaaaaagct
ggaagccgaagtctcgcgcgagattataaacgaaaaagcagaatttctgacctgctggaactg
attgacaacgtccagcgccctgggctgggttaccgtttcgagctcgatccgtgggtgcgctgg
atcgcttcgtttcctccggcggttcgatgcggttaaccaagacttccctgcacggtacggcact
gtcttccgtctgctgcgtcaacacggttttgaggtttctcaggaagcgttcagcggcttcaaa

Figure 37C

Nucleotide sequence of MD09-163 (pET24a-P. alba MEA(+)-TEV)

gacccaaacggcaacttcctggagaacctgaaggaagatatcaaagctatcctgagcctgtacg
agggcagcttcctggctctggaaggcgaaaacatcctggacgaggcgaaaggttttcgcaatctc
tcatctgaaagaactgtctgaagaaaagatcggtaaagagctggcagaacaggtgaaccatgca
ctggaaactgccactgcatcgccgtactcagcgtctggaagcagtatggctctatcgaggcctacc
gtaaaaaggaggacgcgaatcaggttctgctggagctggcaattctggattacaacatgatcca
gtctgtataccagcgtgatctgctgaaacgtcccggtgggtggcgtcgtgtgggtctggcgacc
aaactgcacttttgcctgctgaccgcctgattgagagcttctactgggccgtgggtgtagcattcg
aaccgcaataactccgactgccgtaactccgtcgcaaaaatgttttcttctgtaaccattatcga
cgatatctacgatgtatacggcaccctggacgaactggagctgtttactgatgcagttgagcgt
tgggacgtaaacgccatcaacgacctgccggattacatgaaactgtgctttctggctctgtata
acactattaacgaaatcgctacgacaacctgaaagataaaagtgagaacatcctgccgtatct
gaccaaagcctgggctgacctgtgcaacgctttcctgcaagaagccaagtggctgtacaacaaa
tctactccgacctttgacgactacttcggcaacgcgatggaaatcctcttctggcccgtgcaac
tgggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatcgaaaacctgcaaaa
ataccatgacaccatctctcgtccttcccatatcttccgtctgtgcaatgacctggctagcgcg
tctgcggaaattgcgcgtggtgaaaccgcaaatagcgtttcttgttacatgcgcactaaaggta
tctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacctggaaaaagatgaa
caaggaaaaactgggtggtagcctgttcgcgaaaccgttcgtggaaaccgcgatcaacctggca
cgtcaatctcactgcacttatcataacggcgacgcgcataacctctccggatgagctgacccgca
aacgcgttctgtctgtaatactgaaccgattctgccgtttgaacgcgaaaaacctgtattttca
gggcctcgagcaccaccaccaccactgagatccggctgctaacaaagcccgaaggaagct
gagttggctgctgccaccgctgagcaataactagcataaaccccttggggcctctaaacgggtct
tgaggggttttttgctgaaaggaggaactatatccggat

Figure 38

Amino acid sequence of *P. alba* FL (+) TEV (SEQ ID NO:134)

MRC SVSTENV SFTETETETEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINN
EKA EFLTLLELIDNVQRLGLGYRFESDIRGALDRFVSSGGF DAVTKTSLHGTALSFRLLRQHGF
EVSQEAFSGFKDQNGNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKI
GKELAEQVNHAE LPLHRR TQRLEAVWSIEAYRK KEDANQV LLELA ILDYNMIQSVYQ RDLRET
SRWRRRVGLATKLHFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTI IDDIYDVYGTLD
ELELFTDAVERWDVNAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNA
FLQEAKWLYNKSTPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEE IENLQKYHDTISRPSH
IFRLCNDLASASAEIARGETANSVSCYMR TKGISEELATESVMNLIDETWKKMNKEKLGGSLFA
KPFVETAINLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFERENLYFQGLEHHHHHHH

Figure 39A

Nucleotide sequence of MD09-161 (pET24a-P. alba FL(+)-TEV) (SEQ ID NO:135)

tggcgaatgggacgcgcctgtagcggcgcatataagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccttagcgcgcctcctttcgctttcttcccttctttctcg
ccacggttcgccggtttccccgtcaagctctaaatcgggggctccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaacttgattagggtgatggttcacgtagtgggccatcg
ccctgatagacggtttttcgccctttgacgttggagtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgattataaagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacggtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaataatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaataagg
tatcaagtgagaaatcaccatgagtgcgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatcaacc
aaaccgttattcattcgtgattgcccctgagcgcgagacgaaatacgcgatcgctgttaaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaataacctggaatgctgtttcccggggatcgagtggtgagt
aaccatgcatcatcaggagtagggataaaatgcttgatggtcggaagaggcataaaatccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctaccttgccatgtttcag
aaacaactctggcgcacgggcttcccatacaatcgatagattgtcgacactgattgcccagaca
ttatcgcgagccatttataccatataaaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttcccggtgaatatggctcataaacaccccttgattactgtttatgtaagcaga
cagttttattgttcatgacccaaatcccttaacgtgagttttcgttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgcagataccaaataactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgcctacataacctcgctctgctaatectgttacca
gtggctgctgccagtggcgataaagtcgtgtcttaccgggttggaactcaagacgatagttaccgg
ataaaggcgcagcggtcgggctgaacggggggttcgtgcacacagcccagcttggagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaaagggaga
aaggcggacaggtatccggttaagcggcaggggtcggaacaggagagcgcacaggggagcttccag
ggggaaacgcctggtatctttatagtcctgtcgggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcaggggggaggagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttcttccctgcgttatcccctgattctgtgg
ataaccgtattaccgcctttgagtgagctgataccgctcgccgcagccgaacgaccgagcgcag
cgagtcagtgagcgaaggaagcgaagagcgcctgatgcgggtattttctccttacgcatctgtgc
ggtatttcacaccgcatataggtgcaactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacaccccgccaacac
ccgctgacgcgcctgacgggcttgtctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cgttgagtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaagggcggt

Figure 39B

Nucleotide sequence of MD09-161 (pET24a-P. alba FL(+)-TEV)

tttttcctggttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata
ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaaactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtaggtgttccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcaggcgctgacttccgcgtttccagactttacgaacac
ggaaaccgaagaccattcatgttgttgcacaggtgcagacggtttgcagcagcagtcgcttca
cgttcgctcgcgatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgaccccggtggggcccgccatgccggcgataatggcct
gcttctcgccgaaacgttgggtggcgggaccagtgacgaaggcttgagcgagggcggtgcaagat
tccgaataaccgcaagcgacaggccgatcatcgtcgcgctccagcgaaagcgggtcctcgccgaaa
atgaccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagt
cgggcagcatagtcatgccccgcgcccaccggaaggagctgactgggttgaaggctctcaagg
catcggtcgagatcccggtgcctaatgagtgcgtaacttacattaattgcgttgcgctcactg
ccgcctttccagtcgggaaacctgtcgtgccagctgcattaatgaatcgcccaacgcgcgggga
gaggcggtttgcgtattgggcgccagggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttaccgcctggcctgagagagttgcagcaagcgggtccacgctgggttgcgccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatatccgcaccaacgcgcagcccggactcggtaatggcgcgcattg
cgcccagcgccatctgatcgttggcaaccagcatcgagtggaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atttgattgcgagtgagatatttatgccagccagccagacgcagacgcgcgcgagacagaactta
atgggcccgcctaacagcgcgatttgcgtggtgacccaatgcgaccagatgctccacgccagtcg
cgtaccgtcttcatgggagaaaataactgttgatgggtgtctggtcagagacatcaagaaat
aacgcgggaacattagtgacggcagcttccacagcaatggcatcctgggtcatccagcggatagt
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcgcgtttacaggttcgcac
gccgcttcggttaccatcgacaccaccacgctggcaccagttgatcggcgcgagatttaac
gccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgcgacgcgggttgggaatgtaattcagctccgccatcgccgc
ttccactttttcccgcttttcgcagaaacgtggctggcctgggttaccacgcgggaaacgggtc
tgataagagacaccggcataactctgcgacatcgtataacgttactgggtttcacattcaccaccc
tgaattgactctctccgggcgctatcatgccataccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcgcgcgcaagggaatggtgcatagcaaggagatggcgcccaacagtc
ccggccacggggcctgccaccataccacgcgcgaaacaagcgctcatgagcccgaaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataaggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatccgcgaaattaatcgcac
tactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatgcgttgtagcgtgtccaccgaaaatgtgtctttcaccgaaactgaaa
ccgaagctcgtcgttctgcgaactacgaacctaacagctgggactatgattacctgctgtcctc
cgacacggacgagtcctcgaagtatacaagacaaagcgaaaaagctggaagccgaagttcgt
cgcgagattaataacgaaaaagcagaatttctgacctgctggaactgattgacaacgtccagc
gcctgggcctgggttaccgtttcgagtcgatatccgtgggtgcgctggatcgcttcggttctc
cggcggttcgatgcgtaaccaagacttcctgcacggtagcgcactgtcttccgtctgctg

Figure 39C

Nucleotide sequence of MD09-161 (pET24a-P. alba FL(+)-TEV)

cgtcaacacgggttttgagggtttctcaggaagcgttcagcggccttcaaagacccaaacggcaact
tcttggaagacctgaaggaagatatcaaagctatcctgagcctgtacgaggccagcttcctggc
tctggaaaggcgaaaacatcctggacgaggcgaaaggttttcgcaatctctcatctgaaagaactg
tctgaagaaaagatcggtaaaagagctggcagaacaggtgaaccatgcactggaactgccactgc
atcgccgtactcagcgtctggaagcagtatggtctatcgaggcctaccgtaaaaaggaggacgc
gaatcaggttctgctggagctggcaattctggattacaacatgatccagtctgtataccagcgt
gatctgcgtgaaacgtcccgttgggtggcgtcgtgtgggtctggcgaccaaactgcactttgctc
gtgaccgcctgattgagagcttctactggccgtgggtgtagcattcgaaccgcaatactccga
ctgccgtaactccgtcgcaaaaatgttttcttcgttaaccattatcgacgatatctacgatgta
tacggcacccctggacgaactggagctgtttactgatgcagttgagcgttgggacgtaaacgcc
tcaacgacctgccggttacatgaaactgtgctttctggctctgtataaacactattaacgaaat
cgcctacgacaacctgaaagataaagggtgagaacatcctgccgtatctgaccaaagcctgggt
gacctgtgcaacgctttcctgcaagaagccaaagtggctgtacaacaaatctactccgacctttg
acgactacttcggcaacgcattgaaatcctcttctggcccgctgcaactgggtgttcgcttactt
cgctgtcgtgcagaacattaaaaaggaagagatcgaaaacctgcaaaaataccatgacaccatc
tctcgtccttcccatatcttccgtctgtgcaatgacctggctagcgcgtctgcggaaattgcgc
gtggtgaaaccgcaaatagcgtttcttgttacatgcgcactaaagggtatctccgaagaactggc
taccgaaagcgtgatgaatctgatcgatgaaacctggaaaaagatgaacaaggaaaaactgggt
ggtagcctgttcgcgaaaccgttcgtggaaaccgcgatcaacctggcacgtcaatctcactgca
cttatcataacggcgacgcgcataacctctccggtgagctgacccgcaaacgcgttctgtctgt
aatcaactgaaccgattctgccgtttgaacgcg**gaaaacctgtattttcagggc**ctcgagcaccac
caccaccaccactgagatccggctgctaaacaaagcccgaaggaagctgagttggctgctgcc
ccgtgagcaataactagcataaaccccttggggcctctaaacgggtcttgaggggttttttgc
gaaaggaggaactatatccggat

Figure 40




Figure 41




Figure 42




Figure 43





Figure 44






Figure 45




Figure 46




Figure 47




Figure 48A




Figure 48B




Figure 48C




Figure 49A




Figure 49B




Figure 49C




Figure 50A




Figure 50B




Figure 50C




Figure 51

Amino acid sequence of *P.alba* TRC (-3) (SEQ ID NO:136)

MTEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLELIDNV
QRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFVVSQEAFFSGFKDQNG
NFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVNHAELEP
LHRRTRQLEAVWSIEAYRKEDANQVLELAILDYNMIQSVYQDLRETSRWRRRVGLATKLHF
ARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLDLELEFTDAVERWDVN
AINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNKSTPT
FDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEELENLQKYHDTISRP SHIFRLCNDLASASAEI
ARGETANSVSCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLARQSH
CTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 52A

Nucleotide sequence of pDu47-3 (SEQ ID NO:137)

tggcgaatgggacgcgcctgtagcggcgcat taagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccctagcgcgcgctcctttcgctttcttcccttctttctcg
ccacgttcgcccgtttccccgtcaagctctaaatcgggggctccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaacttgattagggatgatggttcacgtagtgggcatcg
ccctgatagacggtttttcgccctttgacgttggagtccacgttctttaatagtggactcttgt
tccaaactggaacaactcaaccctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaataatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttccctcgtcaaaaaataagg
tatcaagtgagaaatcacatgagtgacgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcacaaacc
aaaccgttattcattcgtgattgcccctgagcgcgagcgaataacgcgatcgtgtttaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatcaacaatattttc
acctgaatcaggatattcttctaataacctggaatgctgttttccggggatcgcagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatcgggttcccatacaatcgatagattgtcgcacctgattgcccagaca
ttatcgcgagccatttataccatataaatcagcatccatgttggaaatttaatcgcggcctag
agcaagacgtttccggttgaatatggctcataaacaccccttgattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcgttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgcccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgcagataccaaataactgtccttctagtgtagccgtagt
aggccaccacttcaagaactctgtagcaccgcctacataacctcgctctgctaactcctgttacca
gtggctgctgccagtggcgataagtcgtgtcttaccgggttggactcaagacgatagttaccgg
ataaggcgcagcggtcgggtgaacggggggttcgtgcacacagcccagcttggagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaagggaga
aaggcggacaggtatccggttaagcggcagggtcggaacaggagagcgcacgagggagcttccag
gggaaacgcctggtatctttatagtcctgtcgggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatccctgattctgtgg
ataaccgtattaccgcctttgagtgagctgataccgctcgcgcgagccgaacgaccgagcgcag
cgagtcagtgcgaggaagcgaagagcgcctgatgcggtattttctccttacgcatctgtgc
ggtatttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacaccgcgaacac
ccgctgacgcgcctgacgggttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cgttgagtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaaggcggt
tttttctgtttggtcactgatgcctccgtgttaagggggatttctgttcatgggggtaatgata

Figure 52B

Nucleotide sequence of pDu47-3

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggttccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcgcaggtcgcagacgttttgcagcagcagtcgcttca
cggttcgctcgcgtatcgggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgcccggcgataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtgcgaaggcttgagcgagggcggtgcaagat
tccgaataccgcaagcgacaggccgcatcatcgtcgcgctccagcgaaagcggctcctcgccgaaa
atgacccagagcgctgcccgcacctgtcctacgagttgcatgataaagaagacagtcataagt
cggcgacgatagtcatgccccgcgcccaccggaaggagctgactgggttgaaggctctcaagg
catcggtcgagatcccggtgctaatgagtgcgtaacttacattaattgcgttgcgctcactg
cccgtttccagtcgggaaacctgtcgtgccagctgcattaatgaatcgcccaacgcgcgggga
gaggcggtttgcgtattgggcgccagggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttcaccgctggccctgagagagttgcagcaagcgggtccacgctgggttgcgccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtctcggtatcg
tcgtatcccactaccgagatatccgcaccaacgcgcagcccggactcggtaatggcgcgcatg
cgcccagcgccatctgatcgttggcaaccagcatcgagtggaacgatgccctcattcagcat
ttgcatgggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atattgattgcgagtgagatatattatgccagccagccagacgcagacgcgcgcgagacagaactta
atgggcccgttaacagcgcgatttgcgtggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaaataactgttgatgggtgtctggtcagagacatcaagaaat
aacgcccgaacattagtgaggcagcttccacagcaatggcatcctggtcatccagcggatagt
taatgatcagcccactgacgcgttgcgcgagaaagattgtgcaccgcccgtttacaggcttcgac
gccgttctggttctaccatcgacaccaccagctggcaccagttgatcggcgcgagatttaatc
gccgcgacaatttgcgacggcgcgtgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgcgcacgcggttgggaatgtaatcagctccgccatcgccgc
ttccactttttcccgcgttttcgcagaaacgtggctggcctggttcaccacgcgggaaacgggtc
tgataagagacaccggcataactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataaccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcccgcgcaagggaatggtgcatgcaaggagatggcgcccaacagtc
ccggccacggggcctgccaccatacccacgcgaaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaatataacgac
tcaactatagggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatgaccgaagctcgtcgttctgcgaactacgaacctaacagctgggact
atgattacctgctgtcctccgacacggacgagtcctcgaagtatacaaagacaaagcgaaaaa
gctggaagccgaagtctcgtcgcgagattaataacgaaaaagcagaatttctgacctgctggaa
ctgattgacaacgtccagcgcctgggcctgggttaccgttccgagctcgatatccgtggtgcgc
tggtatcgcttcgtttcctccggcggcttcgatgcggtaaccaagacttccctgcacggtaacggc
actgtctttccgtctgctgcgtcaacacggttttgaggtttctcaggaagcgttcagcggcttc
aaagacaaaacggcaacttctggagaacctgaagggaagatatcaaagctatcctgagcctgt

Figure 52C

Nucleotide sequence of pDu47-3

acgaggccagcttcctggctctggaaggcgaaaacatcctggacgaggcggaaggttttcgcaat
ctctcatctgaaagaactgtctgaagaaaagatcggtaaagagctggcagaacaggtgaaccat
gcactggaactgccactgcatcgccgtactcagcgtctggaagcagtatggtctatcgaggcct
accgtaaaaaggaggacgcgaatcaggttctgctggagctggcaattctggattacaacatgat
ccagtctgtataccagcgtgatctgctgaaacgtcccgttggaggcgtcgtgtgggtctggcg
accaaactgcactttgctcgtgaccgcctgattgagagcttctactgggccgtgggtgtagcat
tcgaaccgcaatactccgactgccgtaactccgtcgcaaaaatgttttcttctgtaaccattat
cgacgatatctacgatgtatcggcaccctggacgaactggagctgtttactgatgcagttgag
cgttgggacgtaaacgccatcaacgacctgccggattacatgaaactgtgctttctggctctgt
ataacactattaacgaaatcgctacgacaacctgaaagataaagggtgagaacatcctgccgta
tctgaccaaagcctgggtgacctgtgcaacgctttcctgcaagaagccaagtggctgtacaac
aatctactccgacctttgacgactacttcggcaacgcattggaaatcctcttctggcccgtgc
aactggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatcgaaaacctgca
aaaataccatgacaccatctctcgtccttcccatatcttccgtctgtgcaatgacctggctagc
gcgtctgcggaaattgcgcgtggtgaaaccgcaaatagcgtttcttgttacatgcgcaactaaag
gtatctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacctggaaaaagat
gaacaaggaaaaactgggtggtagcctgttcgcgaaaccgttcgtggaaaccgcgatcaacctg
gcacgtcaatctcactgcacttatcataacggcgacgcgcatacctctccggatgagctgaccc
gcaaacgcgttctgtctgtaatcactgaaccgatcttgccgtttgaacgctaaggatccgaatt
cgagctccgtcgacaagcttgccggccgcaactcgagcaccaccaccaccactgagatccggc
tgctaacaaagcccgaaggaagctgagttggctgctgccaccgctgagcaataactagcataa
ccccttggggcctctaaacgggtcttgaggggttttttgctgaaaggaggaactatatccggat

Figure 53

Amino acid sequence of *P.alba* TRC (-4) (SEQ ID NO:138)

METEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLELIDN
VQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFEVSQEA FSGFKDQN
GNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEKIGKELAEQVNHAELEL
PLHRRRTQRLEAVWSIEAYRKKEDANQVLELAILDYNMIQSVYQDLRETSRWRRRVGLATKLH
FARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLDLELELFTDAVERWDV
NAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNKSTP
TFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEEIEENLQKYHDTISRPSHIFRLCNDLASASAE
IARGETANSVSCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLARQS
HCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 54A

Nucleotide sequence of pDu47-4 (SEQ ID NO:139)

tggcgaatgggacgcgcctgtagcggcgcat taagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgcctagcgcgcgtcctttcgctttcttcccttctttctcg
ccacgttcgcgggtttccccgtcaagctctaaatcgggggtccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaaacttgattagggtgatggttcacgtagtgggcatcg
ccctgatagacggtttttcgcccttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaactcaacctatctcggtctattctttgatttataagggttttgcc
gatttcggcctattggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaataatgtatccgctcatgaatttaattcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
ttctgtaatgaaggagaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaataagg
tatcaagtgaagaaatcaccatgagtgcgactgaatccggtgagaatggcaaaagtttatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatcaacc
aaaccgttattcattcgtgattgcgcctgagcgagacgaaatacgcgatcgtgtttaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatcaacaatattttc
acctgaatcaggatattcttctaataacctggaatgctgtttcccggggatcgagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatcgggttcccatacaatcgatagattgtcgacactgattgcccagaca
ttatcgcgagccatttataccatataaaatcagcatccatgttggaaatttaatcgcggcctag
agcaagacgtttcccggtgaatatggctcataacaccccttgattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcggttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggtaactggcttcagcagagcgagataccaaatactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgctacatacctcgtctctgctaactcctgttacca
gtggtgctgccagtgccgataaagtcgtgtcttaccgggttggtactcaagacgatagttaccgg
ataaggcgagcgggtcgggtgaacgggggttcgtgcacacagcccagcttggagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaagggaga
aaggcgagacaggtatccggtgaagcggcaggtcggaacaggagagcgacaggggagcttccag
ggggaacgcctggtatctttatagtctgtcggtttcgccacctctgacttgagcgtcgatt
tttgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatcccctgattctgtgg
ataaccgtattaccgcctttgagtgcgtgataccgctcgccgcagccgaacgaccgagcgag
cgagtcagtgcgaggaagcggaagagcgctgatgcggtattttctccttacgcatctgtgc
ggtattttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacaccgcgaacac
ccgctgacgcgcctgacgggttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcattcacagatgtctgcctgttcatccgcgtccagct
cggttaggtttctccagaagcgttaatgtctggcttctgataaagcgggcatgttaaggcggt
tttttctgtttggtcactgatgcctccgtgaaggggatttctgttcatgggggtaatgata

Figure 54B

Nucleotide sequence of pDu47-4

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacaaactggcggatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggttccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcaggcgctgacttccgcgttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgctcaggtcgcagacgttttgagcagcagtcgcttca
cgttcgctcgcgatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgatacatgcgcacccgtggggccgcatgccggcgataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtgacgaaggcttgagcgagggcggtgcaagat
tccgaataaccgcaagcgacaggccgatacatcgtcgcgctccagcgaaagcggtcctcgccgaaa
atgacccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagt
cggcgacgatagtcatgccccgcgcccaccggaaggagctgactgggttgaaaggctctcaagg
catcggtcgagatcccggtgcctaatgagtgaagtaacttacattaattgcgttgcgctcactg
cccgtttccagtcgggaaacctgtcgtgccagctgcattaatgaatcgccaacgcgcgggga
gaggcggtttgcgtattgggcgcccagggtggtttttctttccaccagtgcagcgggcaacagct
gattgcccttcaccgcctggccctgagagagttgcagcaagcggtccacgctggtttgccccag
caggcgaaaatcctgtttgatggtggttaacggcgggataataacatgagctgtcttcggtatcg
tcgtatccactaccgagataatccgcaccaacgcgcagcccggactcggtaatggcgcgcat
cgccccagcgccatctgatcgttggcaaccagcatcgagtgggaaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atltgatgagagtgagataatgatgccagccagccagacgcagacgcgcgagacagaactta
atgggcccgcataacagcgcgatttgctggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaataatactgttgatgggtgtctggtcagagacatcaagaaat
aacgcgggaacattagtgagggcagcttccacagcaatggcatcctggtcatccagcggatagt
taatgatcagccactgacgcgttgccgcgagaagattgtgcaccgcccgtttacaggcttcgac
gccgttctgttaccatcgacaccaccagctggcaccagttgatcggcgcgagatttaac
gccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgccacgcggttggaatgtaatcagctccgccatcgccgc
ttccacttttcccgcggttttcgcagaaacgtggctggcctggttcaccacgcgggaaacggctc
tgataagagacaccggcataactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttcgggcgctatcatgccataaccgcgaaagggttttgccgcatcagatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggt
gaggccgttgagcaccgcccgcgcaagggaatggtgcatagcaaggagatggcgcccaacagctccc
ccggccacggggcctgccaccatacccacgcggaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgataataggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaatttaatacgac
tcactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatggaaaccgaagctcgtcgttctgcgaactacgaacctaacagctggg
actatgattacctgctgtcctccgacacggacgagtcctatcgaagtatacaaagacaaagcgaa
aaagctggaagccgaagttcgtcgcgagattaataacgaaaaagcagaatttctgacctgctg
gaactgattgacaacgtccagcgccctgggctgggttaccgtttcgagctctgatatccgtgggtg
cgctggatcgcttcgtttcctccggcggttcgatgcggtaaccaagacttccctgcacggtac
ggcactgtctttccgtctgctgcgtcaacacggttttgaggtttctcaggaagcgttcagcggc
ttcaaagacaaaacggcaacttctggagaacctgaaggaagatatcaaagctatcctgagcc

Figure 54C

Nucleotide sequence of pDu47-4

tgtacgaggccagcttcctggctctggaaggcgaaaacatcctggacgaggcgaagggttttcgc
aatctctcatctgaaagaactgtctgaagaaaagatcggtaaagagctggcagaacaggtgaac
catgcaactggaactgccactgcatcgccgtactcagcgtctggaagcagtatggtctatcgagg
cctaccgtaaaaaggaggacgcgaatcaggttctgctggagctggcaattctggattacaacat
gatccagtctgtataccagcgtgatctgctgaaacgtcccgttgggtggcgtcgtgtgggtctg
gcgaccaaactgcactttgctcgtgaccgcctgattgagagcttctactgggccgtgggtgtag
cattcgaaccgcaataactccgactgccgtaactccgtcgcaaaaatgttttcttctcgtaacat
tatcgacgatatctacgatgtatcgccaccctggacgaactggagctgtttactgatgcagtt
gagcgttgggacgtaaacgccatcaacgacctgccggttacatgaaactgtgctttctggctc
tgtataaactattaacgaaatcgctacgacaacctgaaagataaagggtgagaacatcctgcc
gtatctgaccaaagcctgggtgacctgtgcaacgctttcctgcaagaagccaagtggctgtac
aacaatctactccgaccttgacgactacttcggcaacgcatggaaatcctcttctggccgc
tgcaactggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatcgaaaacct
gcaaaaataccatgacaccatctctcgtccttcccatatcttccgtctgtgcaatgacctggct
agcgcgtctgcggaaattgcgcgtggtgaaaccgcaaatagcgtttcttggtacatgcgcacta
aaggatctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacctggaaaaa
gatgaacaaggaaaaactgggtggtagcctgttcgcgaaaccgttcgtggaaaccgcgatcaac
ctggcacgtcaatctcactgcacttatcataacggcgacgcgcataacctctccggatgagctga
cccgcaaacgcgttctgtctgtaatactgaaccgattctgccgtttgaacgctaaggatccga
attcgagctccgtcgacaagcttgcgccgcactcgagcaccaccaccaccaccactgagatcc
ggctgctaacaaagcccgaaaggaaagctgagttggctgctgccaccgctgagcaataactagca
taacccttggggcctctaaacgggtcttgaggggttttttgcgtgaaaggaggaactatatccg
gat

Figure 55

Amino acid sequence of *P.alba* TRC (-5) (SEQ ID NO:140)

MTETEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLELID
NVQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTS LHGTALSFRLLRQHGFVVSQEAFSGFKDQ
NGNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVNHAE
LPLHRRRTQRLEAVWSIEAYRKKEDANQVLELAILDYNMIQSVYQDLRETSRWRRRVGLATKL
HFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTI IDDIYDVYGTLDLELELFTDAVERWD
VNAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNKST
PTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEEENLQKYHDTISRPSHIFRLCNDLASASA
EIARGETANSVSCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLARQ
SHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 56A

Nucleotide sequence of pDu47-5 (SEQ ID NO:141)

tggcgaatgggacgcgcctgtagcggcgcattaagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccttagcgcgcctcctttcgctttcttcccttctttctcg
ccacgttcgcgggtttccccgtcaagctctaaatcgggggctcctttagggttccgatttag
tgctttacggcacctcgacccccaaaaacttgattaggggtgatggttcacgtagtgggcatcg
cctgatagacggtttttcgcccttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaataatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaactcaccgaggcaggtccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttccctcgtcaaaaaataagg
tatcaagtgagaaatcaccatgagtgacgactgaatccggtgagaatggcaaaagttaatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatcaacc
aaaccgttattcattcgtgattgcgcctgagcgcgagaaatacgcgatcgtgtttaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaatacctggaatgctgtttcccggggatcgcagtggtgagt
aaccatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcacgggcttcccatacaatcgatagattgtcgacactgattgcccagaca
ttatcgcgagccatttatacctataaaatcagcatccatgttggaaatttaatecgcggcctag
agcaagacgtttcccggtgaatatggctcataaacaccccttgattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcggtccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgcagataccaaatactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgcctacatacctcgtctgctaatectgttacca
gtggctgctgccagtggcgataagtctgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgcagcggctcgggtgaacggggggttcgtgcacacagcccagcttgagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttccggaagggaga
aaggcggacaggtatccggttaagcggcaggttcggaacaggagagcgcacgagggagcttccag
ggggaacgcctggtatctttatagtctgtcgggtttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatccctgattctgtgg
ataaccgtattaccgcctttgagtgcgtgataccgctcgcgcgcagccgaacgaccgagcgcag
cgagtgcgtgagcgcaggaagcgggaagagcgcctgatgcgggtattttctccttacgcatctgtgc
ggtattttcacaccgcatatatggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacacccgccaaacac
ccgctgacgcgcctgacgggcttgtctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcacgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cgttgagttttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaaggggcgt
tttttctgtttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata

Figure 56B

Nucleotide sequence of pDu47-5

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggtgtccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcaggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcaggtcgcagacggtttgcagcagcagtcgcttca
cggttcgctcgcgtatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgccggcgataatggcct
gcttctcgcgaaacgtttggtggcgggaccagtgacgaaggcttgagcgagggcggtgcaagat
tccgaataccgcaagcgacaggccgatcatcgtcgcgtccagcgaaagcggtcctcgccgaaa
atgaccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagtgc
cggcgacgatagtcatgccccgcgcccaccggaaggagctgactgggttgaaggctctcaaggg
catcggtcgagatcccggtgcctaatagtgagctaaactacattaattgcgttgcgctcactg
cccgtttccagtcgggaaacctgtcgtgccagctgcattaatgaatcgccaaacgcgcgggga
gaggcggtttgcgtattgggcgcgaggggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttcaccgctggcctgagagagttgcagcaagcggtccacgctgggttgcgccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatataccgaccaacgcgcagcccggactcggtaatggcgcgcattg
cgcccagcgccatctgatcgttggcaaccagcatcgagtgggaaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtgga
atttgattgcgagtgagataatttatgccagccagccagacgcgagacgcgcgagacagaactta
atgggcccgcctaacagcgcgattttgctggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaaataatactgttgatgggtgtctggtcagagacatcaagaaat
aacgcgggaacattagtgaggcagcttccacagcaatggcatcctggtcatccagcggtatagt
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggcttcgac
gcccgttctgttctaccatcgacaccaccacgctggcaccagttgatcggcgcgagatttaac
gcccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgtgccacgcggttgggaatgtaatcagctccgccatcgccgc
ttccactttttcccgcgttttcgcagaaacgtggctggcctggttcaccacgcgggaaacgggtc
tgataagagacaccggcataactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataaccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggcggttgagcaccgcccgcgcaagggaatggtgcatgcaaggagatggcgcccaacagtc
ccggccacggggcctgccaccatacccacgcggaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataaggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatacgac
tactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatgactgaaaccgaagctcgtcgttctgcgaactacgaacctaacagct
gggactatgattacctgctgtcctccgacacggacgagtcctatcgaaagtatacaaagacaaagc
gaaaaagctggaagccgaagttcgtcgcgagatttaataacgaaaaagcagaattttctgaccctg
ctggaactgattgacaacgtccagcgctgggcctgggttacggtttcgagctcgatccggtg
gtgcgctggatcgcttcgtttcctccggcggttcgatgcggtaaccaagacttccctgcacgg
tacggcactgtcttccgtctgctgcgtcaacacgggttttgaggtttctcaggaagcggttcagc
ggcttcaaagaccaaaccggaacttccctggagaacctgaaggaagatatcaaagctatcctga

Figure 56C

Nucleotide sequence of pDu47-5

gcctgtacgaggccagcttcctggctctggaaggcgaaaacatcctggacgaggcgaaggtttt
cgcaatctctcatctgaaagaactgtctgaagaaaagatcggtaaagagctggcagaacaggtg
aaccatgcactggaactgccactgcatcgccgtactcagcgtctggaagcagtatggtctatcg
aggcctaccgtaaaaaggaggacgcgaatcaggttctgctggagctggcaattctggattacaa
catgatccagtctgtataccagcgtgatctgcgtgaaacgtcccgttgggtggcgtcgtgtgggt
ctggcgaccaaactgcactttgctcgtgaccgcctgattgagagcttctactgggcccgtgggtg
tagcattcgaaccgcaatactccgactgccgtaaactccgtcgcaaaaatgttttcttctgtaac
cattatcgacgatatactacgatgtatacggcaccctggacgaactggagctgtttactgatgca
gttgagcgttgggacgtaaacgccatcaacgacctgccggattacatgaaactgtgctttctgg
ctctgtataaactattaacgaaatcgccctacgacaacctgaaagataaaggtgagaacatcct
gccgtatctgaccaaagcctgggctgacctgtgcaacgctttcctgcaagaagccaagtggctg
tacaacaaatctactccgacctttgacgactacttcggcaacgcatggaaatcctcttctggcc
cgctgcaactggtgttcgcttacttcgctgtcgtgcagaacattaaaaagggaagagatcgaaaa
cctgcaaaaataccatgacaccatctctcgtccttcccatatcttccgtctgtgcaatgacctg
gctagcgcgtctgcggaaattgcgcgtggtgaaaaccgcaaatagcgtttcttgttacatgcgca
ctaaaggtatctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacctggaa
aaagatgaacaaggaaaaactgggtggtagcctgttcgcgaaaccgttcgtggaaaccgcgac
aacctggcacgtcaatctcactgcacttatcataacggcgacgcgcataacctctccggatgagc
tgaccgcgaaacgcgttctgtctgtaatcactgaaccgattctgccgtttgaacgctaaggatc
cgaattcgagctccgtcgacaagcttgcgcccgactcgagcaccaccaccaccaccactgaga
tccggctgctaacaagcccgaaggaagctgagttggctgctgccaccgctgagcaataacta
gcataaacccttggggcctctaaacgggtcttgaggggttttttgcgtgaaaggaggaactatat
ccggat

Figure 57

Amino acid sequence of *P.alba* TRC (-6) (SEQ ID NO:142)

METETEARRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLELI
DNVQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFEVSQEA FSGFKD
QNGNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVNHAL
ELPLHRRRTQRLEAVWSIEAYRKKEDANQVLELAILDYNMIQSVYQRDLRETSRWRRRVGLATK
LHFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLDELELFTDAVERW
DVNAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNKS
TPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEE IENLQKYHDTISRPSHIFRLCNDLASAS
AEIARGETANSVSCYMRTKGISEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLAR
QSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 58A

Nucleotide sequence of pDu47-6 (SEQ ID NO:143)

tggcgaatgggacgcgcctgtagcggcgcat taagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccttagcgcgcctcctttcgctttcttcccttccctttctcg
ccacgttcgcgggtttccccgtcaagctctaaatcgggggtccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaaacttgattagggtgatggttcacgtagtgggcatcg
ccctgatagacggtttttcgccctttgacgttggagtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttggt
tatttttctaaatacattcaaatatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaataagagaaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttccctcgtcaaaaaataagg
tatcaagtgagaaatcaccatgagtacgactgaatccggtgagaatggcaaaagttaatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatacacc
aaaccgttattcattcgtgattgcgctgagcgagacgaaatacgcgatcgtgtttaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaataacctggaatgctgttttcccggggatcgagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatcgggttcccatacaatcgatagattgtcgacactgattgcccgaca
ttatcgcgagccatttatacccatataaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttcccggtgaatatggctcataacaccccttgattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcggttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgcccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgagataccaaatactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgcctacatacctcgctctgctaatacctgttacca
gtggctgctgccagtggcgataagtcgtgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgagcggtcgggctgaacggggggttcgtgcacacagcccagcttggaagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaagggaga
aaggcggacaggtatccggttaagcggcagggtcggaacaggagagcgacaggggagcttccag
gggaaacgcctggtatctttatagtcctgtcggtttcgccacctctgacttgagcgtcgatt
tttgatgctcgtcagggggcgagcctatggaaaaacgcccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatccctgattctgtgg
ataaccgtattaccgcctttgagtgagctgataccgctcgccgcagccgaacgaccgagcgag
cgagtcagtgagcgaggaagcgggaagagcgccctgatgcggtattttctccttacgcatctgtgc
ggtatttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacacccgccaacac
ccgctgacgcgcctgacgggcttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cggttagtttctccagaagcgttaaatgtctggctctgataaagcgggcatgttaaggcggt
ttttcctggttggtcactgatgcctccgtgtaagggggatttctgttcatgggggtaatgata

Figure 58B

Nucleotide sequence of pDu47-6

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtaggtgttccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcctcaggtcgcagacgttttgcagcagcagtcgcttca
cgttcgctcgcgatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgcccggcgataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtgacgaaggcttgagcgagggcggtgcaagat
tccgaataaccgcaagcgacaggccgcatcatcgctcgcgctccagcgaaagcggtcctcgccgaaa
atgacccagagcgctgcccgcacctgtcctacgagttgcatgataaagaagacagtcataagtg
cggcgacgatagtcatgccccgcgcccaccggaaggagctgactgggttgaaaggctctcaaggg
catcggtcgagatcccggtgcctaatgagtgagctaacttacattaattgcgttgcgctcactg
cccgtttccagtcgggaaacctgtcggtccagctgcattaatgaatcgcccaacgcgcgggga
gaggcggtttgcgtattgggcgcccagggtggtttttcttttccaccagtgcagcgggcaacagct
gattgcccttcaccgcctggcctgagagagttgcagcaagcggtccacgctgggttgccccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatataccgcaccaacgcgcagcccggactcggtaatggcgcgcattg
cgcccagcgccatctgatcggtggcaaccagcatcgagtggaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atattgattgcgagtgagatatttatgccagccagccagacgcagacgcgcgagacagaactta
atgggcccgcctaacagcgcgatttgcgtggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaaataactgttgatgggtgtctggtcagagacatcaagaaat
aacgcgggaacattagtgacggcagcttccacagcaatggcatcctggctcatccagcggatagt
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggttcgcac
gcccgttctgttctaccatcgacaccaccagctggcaccagttgatcggcgcgagatttaate
gccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgcacgcgggttgggaatgtaatcagctccgccatcgccgc
ttccactttttccgcggttttcgcagaaacgtggctggcctgggttcaccacgcgggaaacggtc
tgataagagacaccggcataactctgcgacatcgataaacgttactgggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcccgcgcaagggaatggtgcatgcaaggagatggcgcccaacagtc
ccggccacggggcctgccaccataccacgcggaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataggcgccagcaacgcacctgtggcgccggt
gatgcgggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatacgac
tcactataggggaattgtgagcggataaacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatggaaactgaaaccgaagctcgctcggttctgcgaactacgaacctaa
gctgggactatgattacctgctgtcctccgacacggacgagtcctatcgaagtatacaagacaa
agcgaaaaagctggaagccgaagttcgtcgcgagattaataacgaaaaagcagaatttctgacc
ctgctggaactgattgacaacgtccagcgctgggctgggttaccgtttcgagctctgatatcc
gtggtgcgctggatcgcttcgtttcctccggcggttcgatgcggtaaccaagacttccctgca
cggtagggcactgtcttccgtctgctgcgtcaacacgggttttgaggtttctcaggaagcggtc
agcggcttcaagacccaaacggcaacttctggagaacctgaaggaagatatcaagctatcc

Figure 58C

Nucleotide sequence of pDu47-6

tgagcctgtacgaggccagcttcctggctctggaaggcgaaaacatcctggacgaggcgaaggt
tttcgcaatctctcatctgaaagaactgtctgaagaaaagatcggtaaagagctggcagaacag
gtgaaccatgcaactggaactgccactgcatcgccgtactcagcgtctggaagcagtatggtcta
tcgaggcctaccgtaaaaaggaggacgcgaatcaggttctgctggagctggcaattctggatta
caacatgatccagtctgtataccagcgtgatctgctgaaacgtcccgttggtggcgtcgtgtg
ggtctggcgaccaaactgcactttgctcgtgaccgcctgattgagagcttctactgggccgtgg
gtgtagcattcgaaccgcaatactccgactgccgtaactccgtcgcaaaaatgttttcttctgt
aaccattatcgacgatatctacgatgtatacggcacccctggacgaactggagctgtttactgat
gcagttgagcgttgggacgtaaacgccatcaacgacctgccggattacatgaaactgtgcttct
tggctctgtataaactattaacgaaatcgctacgacaacctgaaagataaaggtgagaacat
cctgccgtatctgaccaaagcctgggctgacctgtgcaacgctttcctgcaagaagccaagtgg
ctgtacaacaaatctactccgacctttgacgactacttcggcaacgcatggaaatcctcttctg
gcccgtgcaactggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatcga
aaacctgcaaaaataccatgacaccatctctcgtccttcccatatcttccgtctgtgcaatgac
ctggctagcgcgtctgcggaaattgcgcgtggtgaaaccgcaaatagcgtttcttgttacatgc
gcactaaaggtatctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacctg
gaaaaagatgaacaaggaaaaactgggtggtagcctgttcgcgaaaccgttcgtggaaaccgcg
atcaacctggcagctcaatctcactgcacttatcataacggcgacgcgcataacctctccgatg
agctgaccgcgcaaacgcgttctgtctgtaatcactgaaccgattctgccgtttgaacgctaagg
atccgaattcgagctccgtcgacaagcttgcggccgcactcgagcaccaccaccaccaccactg
agatccggctgctaacaaagcccgaaggaagctgagttggctgctgccaccgctgagcaataa
ctagcataaacccttggggcctctaaacgggtcttgaggggttttttctgtaaggagggaacta
tatccgat

Figure 59

Amino acid sequence of *P.alba* TRC (-7) (SEQ ID NO:144)

MTETETEARRSANYEPNSWDYDYLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLLEL
IDNVQRLGLGYRFESDIRGALDRFVSSGGFDAVTKTSLHGTALSFRLLRQHGFEVSQEAFSGFK
DQNGNFLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVNHA
LELPLHRRRTQRLEAVWSIEAYRKKEDANQVLELAILDYNMIQSVYQRDLRETSRWRRRVGLAT
KLHFARDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTI IDDIYDVYGTLDELELFTDAVER
WDVNAINDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNK
STPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEEIEENLQKYHDTISRPSHIFRLCNDLASA
SAEIARGETANSVSCYMRTKGISEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLA
RQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 60A

Nucleotide sequence of pDu47-7 (SEQ ID NO:145)

tggcgaatgggacgcgcctgtagcggcgcat taagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccctagcggccgctcctttcgctttcttcccttctttctcg
ccacggttcgcccgtttccccgtcaagctctaaatcgggggtccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaacttgattagggtgatggttcacgtagtgggcatcg
cctgatagacggtttttcgccctttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaat taaacgcgaattttaacaaa
atattaacgtttacaatttcagggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaataatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaataagagaaaaactcaccgaggcaggttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaat tccccctcgtcaaaaaataagg
tatcaagtgaataacacatgagtgacgactgaatccggtgagaatggcaaaagt ttagcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcataaacc
aaaccgttattcattcgtgattgcgcctgagcgcgagcaaatacgcgatcgtgttaaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaatacctggaatgctgttttccgggggatcgagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaat tccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatacgggttcccatacaatcgatagattgtcgcacctgattgcccgaca
ttatcgcgagccatttatacctatataaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttccggttgaatatggctcataaacaccttgtattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcggttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggttaactggcttcagcagagcgcagataccaaatactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgcctacatacctcgctctgctaatacctgttacca
gtggctgctgccagtggcgataagtctgtcttaccgggttggtgactcaagacgatagttaccgg
ataaggcgagcgggtcgggtgaacggggggttcgtgcacacagcccagcttgaggcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaaagggaga
aaggcggacaggtatccggttaagcggcaggggtcggaacaggagagcgcacgagggagcttccag
ggggaacgcctggtatctttatagtcctgtcggtttcgccacctctgacttgagcgtcgatt
tttgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatccctgattctgtgg
ataaccgtattaccgcctttgagtgcgctgataccgctcgccgcagccgaacgaccgagcgcag
cgagtgcgtgagcaggaagcgaagagcgcctgatgcggtattttctccttacgcatctgtgc
ggtatttcacaccgcataatggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacacccgccaacac
ccgctgacgcgcctgacgggttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cgttgagtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaaggcggt
ttttcctgtttggtcactgatgcctccgtgtaagggggtttctgttcatgggggttaataata

Figure 60B

Nucleotide sequence of pDu47-7

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacaactggcggatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtagggttccacagggtagccagcagcatcctgcg
atgcagatccggaacataatgggtgcaggcgctgacttccgcgttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcgcaggtcgcagacgttttgcagcagcagtcgcttca
cggttcgctcgcgtatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgcccggcgataatggcct
gcttctcgcgaaacgtttgggtggcgggaccagtgacgaaggcttgagcgaggcggtgcaagat
tccgaataccgcaagcgacaggccgcatcatcgtcgcgctccagcgaaagcggctcctcgccgaaa
atgacccagagcgctgcccgcacctgtcctacgagttgcatgataaagaagacagtcataagt
cggcgacgatagtcagtcgcccgcgcccaccggaaggagctgactgggttgaaaggctctcaagg
catcggtcgagatcccggtgcctaatgagtgagctaacttacattaattgcgttgcgctcactg
ccgcgttccagtcgggaaacctgtcgtgccagctgcattaatgaatcgccaacgcgcgggga
gaggcggtttgcgtattgggcgccagggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttcaccgcctggcctgagagagttgcagcaagcggccacgctggtttgcccag
caggcgaaaatcctgtttgatgggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatataccgcaccaacgcgcagcccggactcggtaatggcgcgcat
cgcccagcgccatctgatcgttggcaaccagcatcgcagtggaacgatgccctcattcagcat
ttgcatgggttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atttgattgcgagtgagatatttatgccagccagccagacgcagacgcgcgagacagaactta
atgggcccgcctaacagcgcatattgctggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcagggagaaaataatctgttgatgggtgtctggtcagagacatcaagaaat
aacgcccgaacattagtgaggcagcttccacagcaatggcatcctggctcatccagcggatagt
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggcttcgac
gcccgttctgttctaccatcgacaccaccagctggcaccagttgatcggcgcgagatttaac
gcccgcgacaatttgcgacggcgcgtgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgcagttgttgcacgcggttgggaatgtaattcagctccgccatcgccgc
ttccactttttccgcggttttcgagaaacgtggctggcctggttcaccacgcgggaaacggctc
tgataagagacaccggcatactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttcggggcgctatcatgccataaccgcgaaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggt
gaggccgttgagcaccgcccgcgcaaggaatggtgcatagcaaggagatggcgcccaacagtc
ccggccacggggcctgccaccatacccacgcgaaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataggcgccagcaacccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaatataacgac
tcactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatgaccgaaactgaaaccgaagctcgtcgttctgcgaactacgaacct
acagctgggactatgattacctgctgtcctccgacacggacgagtcctatcgaagtatacaaaga
caaagcgaaaaagctggaagccgaagtgcgtcgcgagattaataacgaaaaagcagaatttctg
accctgctggaactgattgacaacgtccagcgctgggcctgggttaccggttctgagctctgata
tccgtggtgcgctggatcgcttcgtttcctccggcggttcgatgcggtaaccaagacttccct
gcacggtacggcactgtctttccgtctgctgcgtcaacacgggttttgaggtttctcaggaagcg
ttcagcggcttcaaagacaaaacggcaacttctggagaacctgaaggaagatatcaaagcta

Figure 60C

Nucleotide sequence of pDu47-7

t c c t g a g c c t g t a c g a g g c c a g c t t c c t g g c t c t g g a a g g c g a a a c a t c c t g g a c g a g g c g a a
g g t t t t c g c a a t c t c t c a t c t g a a a g a a c t g t c t g a a g a a a g a t c g g t a a a g a g c t g g c a g a a
c a g g t g a a c c a t g c a c t g g a a c t g c c a c t g c a t c g c c g t a c t c a g c g t c t g g a a g c a g t a t g g t
c t a t c g a g g c c t a c c g t a a a a g g a g g a c g c g a a t c a g g t t c t g c t g g a g c t g g c a a t t c t g g a
t t a c a a c a t g a t c c a g t c t g t a t a c c a g c g t g a t c t g c g t g a a a c g t c c c g t t g g t g g c g t c g t
g t g g g t c t g g c g a c c a a a c t g c a c t t t g c t c g t g a c c g c c t g a t t g a g a g c t t c t a c t g g g c c g
t g g g t g t a g c a t t c g a a c c g c a a t a c t c c g a c t g c c g t a a c t c c g t c g c a a a a a t g t t t t c t t t
c g t a a c c a t t a t c g a c g a t a t c t a c g a t g t a t a c g g c a c c c t g g a c g a a c t g g a g c t g t t t a c t
g a t g c a g t t g a g c g t t g g g a c g t a a a c g c c a t c a a c g a c c t g c c g g a t t a c a t g a a a c t g t g c t
t t c t g g c t c t g t a t a a c a c t a t t a a c g a a a t c g c c t a c g a c a a c c t g a a a g a t a a a g g t g a g a a
c a t c c t g c c g t a t c t g a c c a a a g c c t g g g c t g a c c t g t g c a a c g c t t t c c t g c a a g a a g c c a a g
t g g c t g t a c a a c a a t c t a c t c c g a c c t t t g a c g a c t a c t t c g g c a a c g c a t g g a a a t c c t c t t
c t g g c c c g c t g c a a c t g g t g t t c g e t t a c t t c g c t g t c g t g c a g a a c a t t a a a a g g a a g a g a t
c g a a a c c t g c a a a a t a c c a t g a c a c c a t c t c t c g t c c t t c c c a t a t c t t c c g t c t g t g c a a t
g a c c t g g c t a g c g c g t c t g c g g a a a t t g c g c g t g g t g a a a c c g c a a a t a g c g t t t c t t g t t a c a
t g c g c a c t a a a g g t a t c t c c g a a g a a c t g g c t a c c g a a a g c g t g a t g a a t c t g a t c g a t g a a a c
c t g g a a a a a g a t g a a c a a g g a a a a c t g g g t g g t a g c c t g t t c g c g a a a c c g t t c g t g g a a a c c
g c g a t c a a c c t g g c a c g t c a a t c t c a c t g c a c t t a t c a t a a c g g c g a c g c g c a t a c c t c t c c g g
a t g a g c t g a c c c g c a a a c g c g t t c t g t c t g t a a t c a c t g a a c c g a t t c t g c c g t t t g a a c g c t a
a g g a t c c g a a t t c g a g c t c c g t c g a c a a g c t t g c g g c c g c a c t c g a g c a c c a c c a c c a c c a c c a
c t g a g a t c c g g c t g c t a a c a a g c c c g a a a g g a a g c t g a g t t g g c t g c t g c c a c c g c t g a g c a a
t a a c t a g c a t a a c c c c t t g g g g c c t c t a a a c g g g t c t t g a g g g g t t t t t t g c t g a a a g g a g g a a
c t a t a t c c g g a t

Figure 61

Amino acid sequence of *P. tremu* TRC (MET) (SEQ ID NO:146)

METRRSANYEPNSWDYDYLLSSDTDESIEVHKDKAKKLEAEVRREINNEKAEFLTLLELIDNVQ
RLGLGYRFESDIRRALDRFVSSGGFDGVTKTSLHGTALSFRLLRQHGFVVSQEA FSGFKDQNGN
FLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVSHALELPL
HRRTQRLEAVWSIEAYRKKEDANQVLELAILDYNMIQSVYQDLRETSRWRRRVGLATKLHFA
RDRLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLDLELFTDAVERWDVNA
INDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNKSTPTF
DDYFGNAWKSSSGPLQLIFAYFAVVQNIKKEEIEENLQKYHDIISRPSHIFRLCNDLASASAEIA
RGETANSVSCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLARQSHC
TYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 62A

Nucleotide sequence of pDu48 (SEQ ID NO:147)

tggcgaatgggacgcgcctgtagcggcgcat taagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccctagcggccgctcctttcgctttcttcccttctttctcg
ccacgttcgcccgtttccccgtcaagctctaaatcgggggctccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaacttgattagggatggttcacgtagtgggcatcg
ccctgatagacggttttcgccccttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaacactcaaccctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctattggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcagggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaataatgtatccgctcatgaattaatcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaaactcaccgaggcaggttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaataagg
tatcaagtgagaaatcaccatgagtgacgaactgaatccggtgagaatggcaaaagttaatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatacacc
aaaccgttattcattcgtgattgcccctgagcgcagacgaaatacgcgatcgtgttaaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatacaaatattttc
acctgaatcaggatattcttctaataacctggaatgctgtttcccggggatcgcagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaattccgtca
gccagtttagtctgaccatctcatctgtaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatacgggttcccatacaatcgatagattgtcgcacctgattgcccgcaca
ttatcgcgagccatttataccatataaatcagcatccatgttggaaatttaatcgcggcctag
agcaagacgtttcccggtgaatatggctcataaacaccccttgtattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcgttccaactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggtaactggcttcagcagagcgcagataccaaataactgtccttctagtgtagecgtagtt
aggccaccacttcaagaactctgtacaccgcctacatacctcgtctctgctaactcctgttacca
gtggctgctgccagtggcgataagtctgtcttaccgggttggtactcaagacgatagttaccgg
ataaggcgcagcgtcgggctgaacggggggttcgtgcacacagcccagcttgaggcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaggggaga
aaggcggacaggtatccggttaagcggcaggggtcggaacaggagagcgcacgaggggagcttccag
gggaaacgcctggtatctttatagtcctgtcggttttcgccacctctgacttgagcgtcgatt
tttgtgatgctcgtcagggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttcttctcctgcgttatcccctgattctgtgg
ataaccgctattaccgcctttgagtgaactgataccgctcgcgcgagccgaacgaccgagcgcag
cgagtcaagtgcgaggaagcggaagagcgccctgatgcggtattttctccttacgcatactgtgc
ggtatttcacaccgcatataggtgcactctcagtacaatctgctctgatgccgcatagtttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacaccgcgaacac
ccgctgacgcgcctgacgggttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcgattcacagatgtctgcctgttcatccgcgtccagct
cggttagtttctccagaagcgttaatgtctggcttctgataaagcgggcatgttaaggcggt
tttttctgtttggtcactgatgcctccgtgtaaggggatttctgttcatgggggttaataata

Figure 62B

Nucleotide sequence of pDu48

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacaactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtaggtgttccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgttgcctcaggtcgcagacgttttgcagcagcagtcgcttca
cgcttcgctcgcgtatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgccatgccggcgataatggcct
gcttctcgcgaaacgtttggtggcgggaccagtgacgaaggcttgagcggaggcggtgcaagat
tccgaataaccgaagcgacaggccgatcatcgtcgcgctccagcgaagcggtcctcgccgaaa
atgaccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagtg
cggcgacgatagtcatgccccgcgcccaccggaaggagctgactgggttgaaggctctcaaggg
catcggtcgagatcccggtgcctaatgagtgcgtaacttacattaattgcgttgcgctcactg
cccgtttccagtcgggaaacctgtcgtgccagctgcattaatgaatcggccaaacgcgcgggga
gaggcggttgcgtattgggcgccagggtggtttttctttccaccagtgcagcgggcaacagct
gattgcccttcaccgcctggcctgagagagttgcagcaagcggtccacgctgggttgcgccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatatccgcaccaacgcgcagcccggactcggtaatggcgcgcattg
cgcccagcgccatctgatcgttggcaaccagcatcgcagtgggaaacgatgccctcattcagcat
ttgcatgggttgttgaaaaccggacatggcaactccagtcgccttcccgttccgctatcggtga
atttgattgcgagtgagatatattatgccagccagccagacgcgagacgcgcgcgagacagaactta
atgggcccgcctaacagcgcgatttgctggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaaataatactgttgatgggtgtctggtcagagacatcaagaaat
aacgccggaacattagtgcaggcagcttccacagcaatggcatcctgggtcatccagcggatagt
taatgatcagcccaactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggcttcgac
gccgttctgttctaccatcgacaccaccacgctggcaccagttgatcggcgcgagatttaac
gccgcgacaatttgcgacggcgcgctgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttgcacgcgggttgggaatgtaattcagctccgccatcgccgc
ttccacttttcccgcgttttcgcagaaacgtggctggcctgggtcaccacgcgggaaacgggtc
tgataagagacaccggcatactctgcgacatcgtataacggttactgggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataaccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcccgcgcaaggaatggtgcatgcaaggagatggcgcccaacagctccc
ccggccacggggcctgccaccatacccacgcggaacaagcgctcatgagcccgaagtggcgag
cccgatcttcccctcggtgatgtcggcgatataaggcgccagcaaccgcacctgtggcgccgggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatcgcac
tactataggggaattgtgagcggataacaattcccctctagaaataattttgtttaactttaa
gaaggagatatacatatgaaacgcgtcgttctgcgaactacgaacctaacagctgggactatg
attacctgctgtcctccgacacggacgagtcctcgaagtacacaaagacaaagcgaaaaagct
ggaagccgaagtctgcgcgagattaataacgaaaaagcagaatttctgacctgctggaactg
attgacaacgtccagcgctgggctgggttaccgtttcgagctcgatataccgtcgtgcgctgg
atcgcttcgtttcctccggcggttcgatggcgtaaccaagacttccctgcacggtacggcact
gtctttccgtctgctgcgtcaacacggttttgagggttctcaggaagcggtcagcgggttcaaa
gacaaaaacggcaacttctggagaacctgaagggaagatatcaaagctatcctgagcctgtacg

Figure 62C

Nucleotide sequence of pDu48

aggccagcttcctggctctggaaggcgaaaacatcctggacgaggcggaagggttttcgcaatctc
tcatctgaaagaactgtctgaagaaaagatcggtaaagagctggcagaacagggtgtcccatgca
ctggaactgccactgcatcgccgtactcagcgtctggaagcagtatggtctatcgaggcctacc
gtaaaaaggaggacgcgaaccagggttctgctggagctggcaattctggattacaacatgatcca
gtctgtataccagcgtgatctgctgaaacgtcccgttgggtggcgtcgtgtgggtctggcgacc
aaactgcactttgctcgtgaccgcctgattgagagcttctactgggccgtgggtgtagcattcg
aaccgcaatactccgactgccgtaactccgtcgcaaaaatgttttcttctgtaaccattatcga
cgatatctacgatgtatacggcaccttgacgaactggagctgttactgatgcagttgagcgt
tgggacgtaaacgccatcaacgacctgccgattacatgaaactgtgcttctggtctctgtata
acactattaacgaaatcgctacgacaacctgaaagataaagggtgagaacatcctgccgtatct
gaccaaagcctgggtgacctgtgcaacgcttctcctgcaagaagccaagtggctgtacaacaaa
tctactccgacctttgacgactacttcggcaacgcatggaaatcctcttctggcccgtgcaac
tgatcttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatcgaaaacctgcaaaa
ataccatgacatcatctctcgtccttcccatatcttccgtctgtgcaatgacctggctagcgcg
tctgcggaaattgcgcgtggtgaaaccgcaaatagcgtttcttggtacatgcgcactaaaggta
tctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacctggaaaaagatgaa
caaggaaaaactgggtggtagcctgttcgcgaaaaccgttcgtggaaaccgcgatcaacctggca
cgtcaatctcactgcacttatcataacggcgacgcgcataacctctccggatgagctgacccgca
aacgcgttctgtctgtaatactgaaccgattctgccgtttgaacgctaaggatccgaattcga
gctccgtcgacaagcttgcggccgcactcgagcaccaccaccaccactgagatccggtgc
taacaaagcccgaaggaagctgagttggctgctgccaccgctgagcaataactagcataaacc
cttggggcctctaaacgggtcttgaggggttttttgcgtgaaaggaggaactatatccggat

Figure 63

Amino acid sequence of *P. tricho* (TRC) (SEQ ID NO:148)

METRRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVRREINNEKAEFLTLLELIDNVQ
RLGLGYRFESDIRRALDRFVSSGGFDAVTKTS LHATALSFRLLRQHGFVVSQEA FSGFKDQNGN
FLENLKEDIKAILSLYEASFLALEGENILDEAKVFAISHLKELSEEKIGKDLAEQVNHAE LPL
HRRTQRLEAVLSIEAYRKKEADQVLELAILDYNMIQSVYQRDLRETSRWRRRVGLATKLHFA
RDLIESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLELELFTNAVERWDVNA
IDDLPDYMKLCFLALYNTINEIAYDNLKEKGENILPYLTAWADLCNAFLQEAKWLYNKSTPTF
DEYFGNAWKSSSGPLQLVFAYFAVQNIKKEE IENLQKYHDIISRPSHIFRLCNDLASASAEIA
RGETANSVSCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAKPFVETAINLARQSHC
TYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 64A

Nucleotide sequence of pDu49 (SEQ ID NO:149)

tggcgaatgggacgcgcctgt agcggcgcat taagcgcggcggtgtggtggttacgcgcagc
gtgaccgctacacttgccagcgccctagcggcgctcctttcgctttcttcccttcttctcg
ccacgttcgcccgtttccccgtcaagctctaaatcgggggctccctttagggttccgatttag
tgctttacggcacctcgacccccaaaaacttgattagggatgatggttcacgtagtgggcatcg
ccctgatagacggtttttcgccctttgacgttgaggtccacgttctttaatagtggactcttgt
tccaaactggaacaactcaaccctatctcggtctattcttttgatttataagggttttgcc
gatttcggcctatttggttaaaaaatgagctgatttaacaaaaatttaacgcgaattttaacaaa
atattaacgtttacaatttcaggtggcacttttcggggaaatgtgcgcggaacccctatttgtt
tatttttctaaatacattcaaatatgtatccgctcatgaattaattcttagaaaaactcatcga
gcatcaaataaaactgcaatttattcatatcaggattatcaataccatatttttgaaaaagccg
tttctgtaatgaaggagaaaactcaccgaggcagttccataggatggcaagatcctggtatcgg
tctgcgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaaataagg
tatcaagtgaagaatcaccatgagtgcgactgaatccggtgagaatggcaaaagtatatgcat
ttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcactcgcatcaacc
aaaccgttattcattcgtgattgcgctgagcgcgagacgaaatacgcgatcgctgttaaaaggac
aattacaaacaggaatcgaatgcaaccggcgaggaacactgccagcgcatcaacaatattttc
acctgaatcaggatattcttctaataacctggaatgctgttttcccggggatcgagtggtgagt
aacatgcatcatcaggagtacggataaaatgcttgatggtcggaagaggcataaatcccgta
gccagtttagtctgaccatctcatctgaacatcattggcaacgctacctttgccatgtttcag
aaacaactctggcgcatcgggcttcccatacaatcgatagattgtcgacactgattgcccgaca
ttatcgcgagccatttataccatataaatcagcatccatgttggaatttaatcgcggcctag
agcaagacgtttcccggtgaatatggctcataacaccccttgattactgtttatgtaagcaga
cagttttattgttcatgacaaaaatcccttaacgtgagttttcggtccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaa
acaaaaaaaccaccgctaccagcgggtggtttgtttgccggatcaagagctaccaactcttttc
cgaaggtaactggcttcagcagagcgcagataccaaatactgtccttctagtgtagccgtagtt
aggccaccacttcaagaactctgtagcaccgcctacatacctcgctctgctaatactgttacca
gtggctgctgccagtggcgataagtctgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgcgagcgtcgggctgaacggggggttcgtgcacacagcccagcttggagcgaacgac
ctacaccgaactgagatacctacagcgtgagctatgagaaaagcggccagcttcccgaagggaga
aaggcggacaggtatccggtgaacggcgagggtcggaacaggagagcgcacaggggagcttccag
ggggaacgcctggtatctttatagtctgtcggttttcgccacctctgacttgagcgtcgatt
tttgatgctcgtcagggggggcgagcctatggaaaaacgccagcaacgcggcctttttacgg
ttcctggccttttgctggccttttgctcacatgttctttcctgcgttatcccctgatctgtgg
ataaccgtattaccgcctttgagtgcgctgataccgctcgccgcagccgaacgaccgagcgcag
cgagtgcgtgagcgcggaagcgcgaagagcgctgatgcggtattttctccttacgcgatctgtgc
ggtatttcacaccgcatatatggtgcactctcagtacaatctgctctgatgccgcatagttaag
ccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccgacaccgcgaacac
ccgctgacgcgcctgacgggcttgctgctcccggcatccgcttacagacaagctgtgaccgt
ctccgggagctgcatgtgtcagaggttttcaccgctcatcaccgaaacgcgcgagggcagctgcgg
taaagctcatcagcgtggtcgtgaagcattcacagatgtctgcctgttcatccgcgtccagct
cgttgagtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaaggggcgt
tttttctgtttggtcactgatgcctccgtgaaggggatttctgttcatgggggtaataata

Figure 64B

Nucleotide sequence of pDu49

ccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcccggttactgg
aacgttgtgagggtaaacactggcggtatggatgcggcgggaccagagaaaaatcactcaggg
tcaatgccagcgcttcgttaatacagatgtaggtgtccacagggtagccagcagcatcctgcg
atgcagatccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacac
ggaaaccgaagaccattcatgttgtgtgctcaggtcgcagacgttttgacgacgagtcgcttca
cgttcgctcgcgatcggtgattcattctgctaaccagtaaggcaaccccgccagcctagccgg
gtcctcaacgacaggagcacgcatcatgcgcacccgtggggccgcatgcccggcgataatggcct
gcttctcgccgaaacgtttggtggcgggaccagtgacgaaggcttgagcgagggcggtgcaagat
tccgaataccgcaagcgacaggccgatcatcgtcgcgctccagcgaaagcggtcctcgccgaaa
atgaccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagtg
cggcgacgatatgcatgccccgcgcccaccggaaggagctgactgggttgaaaggctctcaaggg
catcggtcgagatcccggtgctaatagagtgagctaacttacattaatggcgttgcgctcactg
cccgtttccagtcgggaaacctgtcgtgccagctgcattaatgaatcgccaacgcgcgggga
gaggcggtttgcgtattgggcgcccagggtggtttttcttttaccagtgagacgggcaacagct
gattgcccttaccgctggcctgagagagttgcagcaagcggtccacgctggtttgccccag
caggcgaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccactaccgagatatccgcaccaacgcgcagcccggactcggtaatggcgcgcattg
cgcccagcgccatctgatcgttggaaccagcatcgagtggaacgatgccctcattcagcat
ttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgttccgctatcggtga
atttgattgcgagtgagatatattatgccagccagccagacgcagacgcgcgcgagacagaactta
atgggcccgcataacagcgcgatttgctggtgacccaatgcgaccagatgctccacgcccagtcg
cgtaccgtcttcatgggagaaaataatactgttgatgggtgtctggtcagagacatcaagaaat
aacgcgggaacattagtgaggcagcttccacagcaatggcatcctggtcatccagcggatagtg
taatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcccgtttacaggcttcgac
gccgcttcgttctaccatcgacaccaccacgctggcaccagttgatcggcgcgagatttaate
gccgcgacaatttgcgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacg
actgtttgcccgccagttgttggtgcccacgcgggttgggaatgtaattcagctccgccatcgccgc
ttccactttttcccgcgttttcgcagaaacgtggctggcctggttcaccacgcgggaaacgggtc
tgataagagacaccggcatactctgcgacatcgtataacgttactggtttcacattcaccaccc
tgaattgactctcttccgggcgctatcatgccataccgcgaaagggttttgcgccattcgatggt
gtccgggatctcgacgctctcccttatgcgactcctgcattaggaagcagcccagtagtaggtt
gaggccgttgagcaccgcccggcaaggaatggtgcatgcaaggagatggcgcccaacagtcctc
ccggccacggggcctgccaccataccacgcgcgaaacaagcgctcatgagcccgaagtggcgag
cccgatcttccccatcggtgatgtcggcgatataaggcgccagcaaccgcacctgtggcgccggt
gatgccggccacgatgcgtccggcgtagaggatcgagatctcgatcccgcgaaattaatacgac
tcactataggggaattgtgagcggataacaatccccctctagaaataattttgtttaactttaa
gaaggagatatacatatgaaacgcgtcgttctgcgaactacgaacctaacagctgggactatg
attacctgctgtcctccgacacggacgagtcctatcgaagtatacaaagacaaagcgaaaaagct
ggaagccgaagtctcgtcgcgagattaataacgaaaaagcagaatttctgacctgctggaactg
attgacaacgtccagcgctgggcctgggttaccgtttcgagtcgatatccgtcgtgcgctgg
atcgcttcgttccctccggcggttcgatgcggttaaccaagacttcctgcacgcgacggcact
gtctttccgtctgctgcgtcaacacgggttttgaggtttctcaggaagcggttcagcggcttcaa
gacaaaaacggcaacttctggagaacctgaaggaagatatcaaagctatcctgagcctgtacg

Figure 64C

Nucleotide sequence of pDu49

aggccagcttcctggctctggaaggcgaaaacatcctggacgaggcgaagggttttcgcaatctc
tcatctgaaagaactgtctgaagaaaagatcggtaaagatctggcagaacagggtgaaccatgca
ctggaactgccactgcatcgccgtactcagcgtctggaagcagtactgtctatcgaggcctacc
gtaaaaaggaggacgcggatcaggttctgctggagctggcaattctggattacaacatgatcca
gtctgtataccagcgtgatctgcgtgaaacgtcccgttggaggcgtcgtgtgggtctggcgacc
aaactgcactttgctcgtgaccgcctgattgagagcttctactgggccgtgggtgtagcattcg
aaccgcaataactccgactgccgtaactccgtcgcaaaaatgttttcttctgtaaccattatcga
cgatatctacgatgtatacggcaccctggacgaactggagctgtttactaacgcagttgagcgt
tgggacgtaaacgccatcgacgatctgccggattacatgaaactgtgcttctggtctgtata
acactattaacgaaatcgctacgacaacctgaaagaaaaagggtgagaacatcctgccgtatct
gaccaaagcctgggctgacctgtgcaacgctttcctgcaagaagccaagtggctgtacaacaaa
tctactccgacctttgacgaatacttcggcaacgcattgaaatcctcttctggcccgctgcaac
tgggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatcgaaaacctgcaaaa
ataccatgacatcatctctcgtccttcccatatcttccgtctgtgcaatgacctggctagcgcg
tctgcggaaattgcgcgtggtgaaaccgcaaatagcgtttcttgttacatgcgcactaaaggta
tctccgaagaactggctaccgaaagcgtgatgaatctgatcgatgaaacctggaaaaagatgaa
caaggaaaaactgggtggtagcctgttcgcgaaaccgttcgtggaaaccgcgatcaacctggca
cgtcaatctcactgcacttatcataacggcgacgcgcataacctctccggatgagctgacctgca
aacgcgttctgtctgtaatactgaaccgattctgccgtttgaacgctaaggatccgaattcga
gctccgtcgacaagcttgcggccgcactcgagcaccaccaccaccactgagatccggctgc
taacaaagcccgaaggaagctgagttggctgctgccaccgctgagcaataactagcataaacc
cttggggcctctaaacgggtcttgagggggttttttgcgtgaaaggaggaactatatccggat

Figure 65

Amino acid sequence of kudzu TRC (MEA) (SEQ ID NO:150)

MEARRSANYQPNLWNFEFLQSLLENDLKVEKLEEKATKLEEEVRCMINRVDTQPLSLLELIDDVQ
RLGLTYKFEKDI IKALENIVLLDENKKNKSDLHATLSFRLLRQHGFEVSQDVFERFKDKEGGF
SGELKGDVQGLLSLYEASYLGFEGENLLEEARTFSITHLKNNLKEGINTKVAEQVSHALELPYH
QRLHRLEARWFLDKYEPKEPHHQLLELAKLDFNMVQTLHQKELQDLSRWWTEMGLASKLDFVR
DRLMEVYFWALGMAPDPQFGECKAVTKMFGLVTIIDDVYDVYGTLDLQLFTDAVERWDVNAI
NTLPDYMKLCFLALYNTVNDTSYSILKEKGHNLSYLTCSWRELCKAFLQEAKWSNNKIIPAFS
KYLENASVSSSGVALLAPSYFSVCQQQEDISDHALRSLTDFHGLVRSSCVIFRLCNDLATSAAE
LERGETTNSII SYMHENDGTSEEQAREELRKLIDAEWKKMNRERVSDSTLLPKAFMEI AVNMAR
VSHCTYQYGDGLGRPDYATENRIKLLLLIDPFPINQLMYV

Figure 66A

Nucleotide sequence of pDu50 (SEQ ID NO:151)

tcgagcaccaccaccaccaccactgagatccggctgctaacaaagcccgaaggaagctgagtt
ggctgctgccaccgctgagcaataactagcataaacccttggggcctctaaacgggtcttgagg
ggttttttgctgaaaggagggaactatatccggattggcgaatgggacgcgcctctagcggcgc
attaagcgcggcgggtgtggtggttacgcgcagcgtgaccgctacacttgccagcgccttagcg
cccgtctctttcgttttcttcccttcttctcgccacgttcgccggtttccccgtcaagctc
taaatacggggggtccctttagggttccgatttagtgctttacggcacctcgaccccaaaaaact
tgattaggggtgatggttcacgtagtgggcatcgccctgatagacggtttttcgccctttgacg
ttggagtccacgttctttaatagtggactcttgttccaaactggaacaacactcaaccctatct
cggctctattcttttgatttataagggattttgccgatttcggcctattggttaaaaaatgagct
gatttaacaaaaatttaacgcgaattttaacaaaatattaacgtttacaatttcaggtggcact
tttcggggaaatgtgcgcggaacccctatttgtttatttttctaaatacattcaaatatgtatc
cgctcatgaattaattcttagaaaaactcatcgagcatcaaatagaaactgcaatttattcatat
caggattatcaataccatatttttgaaaaagccgtttctgtaatgaaggagaaaaactcaccgag
gcagttccataggtatggcaagatcctggtatcggctctgcgattccgactcgtccaacatcaata
caacctattaatttcccctcgtcaaaaataaggttatcaagtgagaaatcaccatgagtgacga
ctgaatccggtgagaatggcaaaagtttatgcatttcttccagacttgttcaacaggccagcc
attacgctcgtcatcaaaatcactcgcatcaaccaaaaccgttatcattcgtgattgcgcctga
gcgagacgaaatacgcgatcgctgttaaaaggacaattacaaacaggaatcgaatgcaaccggc
gcaggaacactgccagcgcatacaaatattttcacctgaatcaggatattcttctaataacctg
gaatgctgttttcccggggatcgagtggtgagtaacatgcatacatcaggagtacggataaaa
tgcttgatggtcggaagaggcataaattccgtcagccagtttagtctgacctctcatctgtaa
catcattggcaacgctacctttgccatgtttcagaaacaactctggcgcatcgggcttcccata
caatcgatagattgtgcacactgattgcccgcattatcgcgagcccatttatacccatataaa
tcagcatccatgttggaatttaacgcggcctagagcaagacgtttcccgttgaaataggtca
taacaccccttgattactgtttatgtaagcagacagttttattgttcatgacaaaaatccctt
aacgtgagttttcgttccactgagcgtcagaccccgtagaaaaagatcaaaggatcttcttgaga
tcctttttttctgcgcgtaatctgctgcttgcaaacaaaaaaaccaccgctaccagcgggtggtt
tgtttgccggatcaagagctaccaactcttttccgaaggtaactggcttcagcagagcgcaga
taccaaatactgtccttctagtgtagccgtagttaggccaccacttcaagaactctgtagcacc
gcctacataacctcgctctgctaatcctgttaccagtggtgctgccagtgccgataagtcgtgt
cttaccgggttggaactcaagacgatagttaccggataaaggcgcagcggtcgggctgaacggggg
gttcgtgcacacagcccagcttgagcgcaacgcctacaccgaactgagataacctacagcgtga
gctatgagaaaagcgccacgcttcccgaaggagaaaaggcggacaggtatccggtgaacgggcagg
gtcggaaacaggagagcgcacgagggagcttccagggggaaacgcctggtatctttatagtcctg
tcgggtttcgccacctctgacttgagcgtcgatttttgtgatgctcgtcaggggggaggagcct
atggaaaaacgccagcaacgcggcctttttacggttcctggccttttgctggccttttgctcac
atgttctttcctgcgttatccctgattctgtggataaacgtattaccgcctttgagtgagctg
ataccgctcgccgcagccgaacgaccgagcgcagcagtcagtgagcaggaagcgggaagagcg
cctgatgcgggtattttctccttacgcactgtgcgggtatttcacaccgcataataggtgcactc
tcagtacaatctgctctgatgccgcatagttaagccagtatatactccgctatcgctacgtgac
tgggtcatggctgcgccccgacacccgcaacacccgctgacgcgccttgacggggttgtctgc
tcccggcatccgcttacagacaagctgtgaccgtctccgggagctgcatgtgtcagaggttttc
accgtcatcaccgaaacgcgcgaggcagctgcggtaaaagctcatcagcgtggtcgtgaagcgat

Figure 66B

Nucleotide sequence of pDu50

tccacagatgtctgcctgttcatccgcgtccagctcgttgagtttctccagaagcgttaatgtct
ggcttctgataaagcgggcatgttaagggcggtttttcctgttgggtcactgatgcctccgt
gtaagggggatttctgttcatgggggtaatgataccgatgaaacgagagaggatgctcacgata
cgggttactgatgatgaacatgcccggttactggaacgttgtgagggtaaacaactggcggtat
ggatgcggcgggaccagagaaaaatcactcaggggtcaatgccagcgcttcgttaatacagatgt
aggtgttccacagggtagccagcagcatcctgcgatgcagatccggaacataatggtgcagggc
gctgacttccgcgtttccagactttacgaaacacggaaaccgaagaccattcatgttgttgc
aggtcgcagacgttttgcagcagcagtcgcttcacgttcgctcgcgtatcggtgattcattctg
ctaaccagtaaggcaaccccgccagcctagccgggtcctcaacgacaggagcacgatcatgcgc
accggtggggccgcatgcccgcgataatggcctgcttctcgccgaaacgtttggggcgggac
cagtgacgaaggcttgagcgaggcggtgcaagattccgaataccgcaagcgacaggccgatcat
cgtcgcgctccagcgaaagcgggtcctcgccgaaaatgacccagagcgctgccggcacctgtcct
acgagttgcatgataaagaagacagtcataagtgcggcgacgatagtcatgccccgcgcccacc
ggaaggagctgactgggttgaaggctctcaaggcgatcggtcgagatcccgggtgcctaatgagt
gagctaaacttacattaattgcgttgcgctcactgcccgctttccagtcgggaaacctgtcgtgc
cagctgcattaatgaatcggccaacgcgcggggagaggcggtttgcgtattgggcgccagggtg
gtttttcttttaccagtgagacgggcaacagctgattgcccttcaccgcctggccctgagaga
gttgcagcaagcgggtccacgctggtttgcccagcaggcgaaaaatcctgtttgatgggtggttaa
cggcgggatataacatgagctgtcttcgggtatcgtcgtatcccactaccgagatataccgacca
acgcgcagcccggactcggtaatggcgcgcattgcccccagcgccatctgatcgttggcaacca
gcategcagtggaacgatgccctcattcagcatttgcatgggttgttgaaaaccggacatggc
actccagtcgccttcccgttccgctatcggctgaatttgattgcgagtgagatatttatgccag
ccagccagacgcagacgcgcgagacagaacttaatgggcccgttaacagcgcgatttgctggt
gacccaatgcgaccagatgctccacgcccagtcgcgtaccgtcttcatgggagaaaaataact
gttgatgggtgtctggtcagagacatcaagaaataacgcgggaacattagtgaggcagcttcc
acagcaatggcatcctggtcatccagcggatagttaatgatcagcccactgacgcgttgcgcga
gaagattgtgcaccgcgcgtttacaggcttcgacgcgcgttcgttctaccatcgacaccaccac
gctggcacccagttgatcggcgcgagatttaatcgccgcgacaatttgcgacggcgcggtgcagg
gccagactggaggtggcaacgccaatcagcaacgactgtttgcccgcagttgttgtgccacgc
ggttgggaatgtaatcagctccgccatcgccgcttccacttttcccgcgttttcgcagaaac
gtggctggcctggttaccacgcgggaaacggtctgataagagacaccggcatactctgcgaca
tcgtataacgttactggtttcacattcaccacctgaattgactctcttccgggcgctatcatg
ccataccgcgaaagggtttgcgccattcgatggtgtccgggatctcgacgctctcccttatgcg
actcctgcattaggaagcagcccagtagtaggttgaggccgttgagcaccgcccgcgcaaggaa
tggtgcatgcaaggagatggcgcccaacagtcctcccgccacggggcctgccaccatacccacg
ccgaaacaagcgtcatgagcccgaagtggcgagcccgatcttccccatcgggtgatgtcggcga
tataggcgccagcaaccgcacctgtggcgccggtgatgccggccacgatgcgtccggcgtagag
gatcgagatctcgatcccgcgaaattaatacgaactcactataggggaattgtgagcggataaca
attcccctctagaaataattttgtttaactttaagaaggagatataccatggaagctcgtcgtt
ccgcaaaactatcagccaaacctgtggaatttcgaattcctgcaatccctggagaacgacctgaa
agtggaaaagctggaggagaaagcgaccaaactggaggaagaagttcgtgcgatgatcaaccgt
gtagacaccagccgctgtccctgctggagctgatcgacgatgtgcagcgccctgggtctgacct
acaaatttgaaaaagacatcattaagccctggaaaacatcgtactgctggacgaaaacaaaaa

Figure 66C

Nucleotide sequence of pDu50

gaacaaatctgacctgcacgcaaccgctctgtctttccgtctgctgcgtcagcacggtttcgag
gtttctcaggatgtttttgagcggttcaaggataaagaagggtggtttcagcggggaactgaaag
gtgacgtccaaggcctgctgagcctgtatgaagcgtcttacctgggtttcgagggtgagaacct
gctggaggaggcgcgtaacctttccatcacccacctgaagaacaacctgaaagaaggcattaat
accaaggttgagaacaagtgaaccacgacctggaactgccatataccagcgtctgcaccgtc
tggaggcacgttgggttccctggataaatacgaaccgaaagaaccgcatcaccagctgctgctgga
gctggcgaagctggattttaacatggtacagacctgcaccagaaagagctgcaagatctgtcc
cgctggtggaccgagatggcctggctagcaaacctggattttgtacgcgaccgacctgatggaag
tttattttctgggcactgggtatggcgccagaccgcagtttggatgaatgtcgaaaagctgttac
taaaatgtttggtctggtgacgatcatcgatgacgtgtatgacgtttatggcactctggacgaa
ctgcaactgttcaccgatgctgtagagcgtgggacgttaacgctattaacacctgccggact
atatgaaactgtgtttccctggcactgtacaacaccttaacgacacgtcctattctattctgaa
agagaaaggtcataacaacctgtcctatctgacgaaaagctggcgtgaactgtgcaaagccttt
ctgcaagaggcgaaatggtccaacaacaaaattatcccggctttctccaagtacctggaaaacg
ccagcgtttcctcctccgggtgtagcgtgctggcgccgtcttactttccgtatgccagcagca
ggaagacatctccgaccacgcgctgcgttcctgaccgacttccatggtctggtgcgttctagc
tgcgttatcttcgcctgtgcaacgatctggccacctctgcggcgagctggaacgtggcgaga
ctaccaattctatcattagctacatgcacgaaaacgatggtaccagcgaggaacaggcccgga
agaactgcgtaaactgatcgacgccgaatggaaaaagatgaatcgtgaacgcgttagcgactcc
acctgctgcctaaagcgttcatggaaatcgagtttaacatggcacgtgtttcccactgcacct
accagtatggcgatggtctgggtcgcccagactacgcgactgaaaaccgcatcaaactgctgct
gattgacctttcccgattaacagctgatgtatgtctaactgcagggatccgaattcgagctc
cgtcgacaagcttgcgccgcac

Figure 67

Amino acid sequence of kudzu TRC (-4) (SEQ ID NO:152)

MEHNSRRSANYQPNLWNFEFLQSLLENDLKVEKLEEKATKLEEEVRCMINRVDTQPLSLLELIDD
VQRLGLTYKFEKDI IKALENIVLLDENKKNKSDLHATLSFRLLRQHGFVSQDVFERFKDKEG
GFSGELKGDVQGLLSLYEASYLGFEGENLLEEARTFSITHLKNNLKEGINTKVAEQVSHALELP
YHQLRHRLEARWFLDKYEPKEPHHQLLLELAKLDFNMVQTLHQKELQDLSRWWTGLASKLDF
VRDRLEMEVYFWALGMAPDPQFGECRKAVTKMFGGLVTIIDDVYDVYGTLDLQLFDAVERWDVN
AINTLPDYMKLCFLALYNTVNDTSYSILKEKGHNLSYLTKSWRELCKAFLQEAKWSNNKIIPA
FSKYLENASVSSSGVALLAPSYFSVCQQQEDISDHALRSLTDFHGLVRSSCVIFRLCNDLATSA
AELERGETTNSIISYMHENDGTSEEQAREELRKLIDAEWKKMNRERVSDSTLLPKAFMEIAVNM
ARVSHCTYQYGDGLGRPDYATENRIKLLLIDPFPINQLMYV

Figure 68A

Nucleotide sequence of pDu50-4 (SEQ ID NO:153)

t c g a g c a c c a c c a c c a c c a c c a c t g a g a t c c g g c t g c t a a c a a g c c c g a a g g a a g c t g a g t t
g g c t g c t g c c a c c g c t g a g c a a t a a c t a g c a t a a c c c c t t g g g g c c t c t a a a c g g g t c t t g a g g
g g t t t t t g c t g a a a g g a g g a a c t a t a t c c g g a t t g g c g a a t g g g a c g c g c c c t g t a g c g g c g c
a t t a a g c g c g g c g g g t g t g g t g g t t a c g c g c a g c g t g a c c g c t a c a c t t g c c a g c g c c c t a g c g
c c c g c t c c t t t c g c t t t c t t c c c t t c c t t t c t c g c c a c g t t c g c c g g c t t t c c c c g t c a a g t c
t a a a t c g g g g g c t c c c t t a g g g t t c c g a t t t a g t g c t t t a c g g c a c c t c g a c c c c a a a a a a c t
t g a t t a g g g t g a t g g t t c a c g t a g t g g g c c a t c g c c c t g a t a g a c g g t t t t t c g c c c t t t g a c g
t t g g a g t c c a c g t t c t t t a a t a g t g g a c t c t t g t t c c a a a c t g g a a c a c a c t c a a c c c t a t c t
c g g t c t a t t c t t t g a t t t a a g g g a t t t t g c c g a t t t c g g c c t a t t g g t t a a a a a t g a g c t
g a t t t a a c a a a a t t t a a c g c g a a t t t t a a c a a a t a t t a a c g t t t a c a a t t t c a g g t g g c a c t
t t t c g g g g a a t g t g c g c g g a a c c c c t a t t t g t t t a t t t t c t a a a t a c a t t c a a a t a t g t a t c
c g t c a t g a a t t a a t t c t t a g a a a a c t c a t c g a g c a t c a a a t g a a a c t g c a a t t t a t t c a t a t
c a g g a t t a t c a a t a c c a t a t t t t t g a a a a g c c g t t t c t g t a a t g a a g g a g a a a a c t c a c c g a g
g c a g t t c c a t a g g a t g g c a a g a t c c t g g t a t c g g t c t g c g a t t c c g a c t c g t c c a a c a t c a a t a
c a a c c t a t t a a t t t c c c c t c g t c a a a a a t a a g g t t a t c a a g t g a g a a t c a c c a t g a g t g a c g a
c t g a a t c c g g t g a g a a t g g c a a a g t t t a t g c a t t t c t t t c c a g a c t t g t t c a a c a g g c c a g c c
a t t a c g t c g t c a t c a a a a t c a c t c g c a t c a a c c a a a c c g t t a t t c a t t c g t g a t t g c g c c t g a
g c g a g a c g a a a t a c g c g a t c g c t g t t a a a a g g a c a a t t a c a a c a g g a a t c g a a t g c a a c c g g c
g c a g g a a c a c t g c c a g c g c a t c a a c a a t a t t t t c a c c t g a a t c a g g a t a t t c t t c t a a t a c c t g
g a a t g c t g t t t t c c c g g g g a t c g c a g t g g t g a g t a a c c a t g c a t c a t c a g g a g t a c g g a t a a a a
t g c t t g a t g g t c g g a a g a g c a t a a a t t c c g t c a g c c a g t t t a g t c t g a c c a t c t c a t c t g t a a
c a t c a t t g g c a a c g t a c c t t t g c c a t g t t t c a g a a c a a c t c t g g c g c a t c g g g c t t c c c a t a
c a a t c g a t a g a t t g t c g c a c c t g a t t g c c c g a c a t t a t c g c g a g c c c a t t t a t a c c c a t a t a a a
t c a g c a t c c a t g t t g g a a t t a a t c g c g g c c t a g a g c a a g a c g t t t c c c g t t g a a t a t g g c t c a
t a a c a c c c c t t g t a t t a c t g t t t a t g t a a g c a g a c a g t t t t a t t g t t c a t g a c c a a a a t c c c t t
a a c g t g a g t t t t c g t t c c a c t g a g c g t c a g a c c c c g t a g a a a g a t c a a a g g a t c t t c t t g a g a
t c c t t t t t t t c t g c g c g t a a t c t g c t g c t t g c a a c a a a a a a a c c a c c g c t a c c a g c g g t g g t t
t g t t t g c c g g a t c a a g a g c t a c c a a c t c t t t t t c c g a a g g t a a c t g g c t t c a g c a g a g c g c a g a
t a c c a a a t a c t g t c c t t c t a g t g t a g c c g t a g t t a g g c c a c c a c t t c a a g a a c t c t g t a g c a c c
g c c t a c a t a c c t c g c t c t g c t a a t c c t g t t a c c a g t g g c t g c t g c c a g t g g c g a t a a g t c g t g t
c t t a c c g g g t t g g a c t c a a g a c g a t a g t t a c c g g a t a a g g c g c a g c g g t c g g g c t g a a c g g g g g
g t t c g t g c a c a g c c c a g c t t g g a g c g a a c g a c c t a c a c c g a a c t g a g a t a c c t a c a g c g t g a
g c t a t g a g a a g c g c c a c g c t t c c c g a a g g g a g a a a g g c g g a c a g g t a t c c g g t a a g c g g c a g g
g t c g g a a c a g g a g a g c g c a c g a g g g a g c t t c c a g g g g g a a a c g c c t g g t a t c t t t a t a g t c c t g
t c g g g t t t c g c c a c c t c t g a c t t g a g c g t c g a t t t t t g t g a t g c t c g t c a g g g g g c g g a g c c t
a t g g a a a a c g c c a g c a a c g c g g c c t t t t a c g g t t c c t g g c c t t t t g c t g g c c t t t t g c t c a c
a t g t t c t t t c c t g c g t t a t c c c c t g a t t c t g t g g a t a a c c g t a t t a c c g c c t t t g a g t g a g c t g
a t a c c g c t c g c c g a g c c g a a c g a c c g a g c g c a g c g a g t c a g t g a g c g a g g a a g c g g a a g a g c g
c c t g a t g c g g t a t t t t c t c c t a c g c a t c t g t g c g g t a t t t c a c a c c g c a t a t a t g g t g c a c t c
t c a g t a c a a t c t g c t c t g a t g c c g c a t a g t t a a g c c a g t a t a c a c t c c g c t a t c g t a c g t g a c
t g g g t c a t g g c t g c g c c c g a c a c c c g c c a a c a c c c g c t g a c g c g c c c t g a c g g g c t t g t c t g c
t c c g g c a t c c g c t t a c a g a c a a g c t g t g a c c g t c t c c g g g a g c t g c a t g t g t c a g a g g t t t t c
a c c g t c a t c a c c g a a a c g c g c g a g g c a g c t g c g g t a a a g c t c a t c a g c g t g g t c g t g a a g c g a t

Figure 68B

Nucleotide sequence of pDu50-4

tacagatgtctgcctgttcacccgctccagctcgttgagtttctccagaagcgttaatgtct
ggcttctgataaagcgggcatgttaagggcgggtttttcctgtttggcactgatgcctccgt
gtaagggggatttctgttcacggggtaatgataccgatgaaacgagagaggatgctcacgata
cgggttactgatgatgaacatgcccggttactggaacgttgtgagggtaaacactggcggat
ggatgcggcgggaccagagaaaaatcactcagggccaatgccagcgcttcgttaatacagatgt
aggtgttccacagggtagccagcagcatcctgcgatgcagatccggaacataatggtgcagggc
gctgacttccgcgtttccagactttacgaaacacggaaccgaagaccattcatgttgctc
aggtcgagacgttttgagcagcagtcgcttcacgttcgctcgctatcggtgattcattctg
ctaaccagtaaggcaaccccgccagcctagccgggtcctcaacgacagagacagcatcatgcgc
acccgtggggccgcatgccggcgataatggcctgcttctcgccgaaacgtttggtggcgggac
cagtgacgaaggcttgagcagggcggtgcaagattccgaataccgcaagcgacagccgatcat
cgtcgcgctccagcgaaagcgggtcctcgccgaaaatgacccagagcgctgccggcacctgtcct
acgagttgcatgataaagaagacagtcataagtgcggcgacgatagtcatgccccgcgcccacc
ggaaggagctgactgggtgaaggctctcaaggcatcggtcgagatcccgggtgcctaatgagt
gagctaacttacattaattgcgttgcgctcactgcccgtttccagtcgggaaacctgtcgtgc
cagctgcattaatgaatcggccaacgcgcggggagagcggtttgcgtattgggcgccagggtg
gtttttctttcaccagtgagacgggcaacagctgattgcccttcaccgcctggccctgagaga
gttgagcaagcgggtccacgctggtttgccccagcaggcgaaaatcctgtttgatggtggttaa
cggcgggatataacatgagctgtcttcgggtatcgtcgtatcccactaccgagatatccgcacca
acgcgcagcccgactcggtaatggcgcgcatgccccagcgccatctgatcgttggaacca
gcatcgagtggaacgatgccctcattcagcatttgcatggtttgttgaaaaccggacatggc
actccagtcgccttcccgttccgctatcggtgaatttgattgagagtgagatatatatgccag
ccagccagacgcagacgcgcgagacagaactaatgggcccgttaacagcgcgatttgctggt
gacccaatgcgaccagatgctccacgcccagtcgctaccgtcttcatgggagaaaaataact
gttgatgggtgtctggtcagagacatcaagaaataacgccggaacattagtgcaggcagcttcc
acagcaatggcatcctggtcatccagcggatagttaatgatcagcccactgacgcgttgcgcg
gaagattgtgacccgcgctttacaggcttcgacgccgttcgttctaccatcgacaccaccac
gctggcaccagttgatcggcgcgagatttaatcgccgcgacaatttgcgacggcgcggtgcagg
gccagactggaggtggcaacgccaatcagcaacgactgtttgcccgccagttgttggtgccacgc
ggttggaatgtaattcagctccgccatcgccgcttccacttttcccgcttttcgcagaaac
gtggctggcctggttcaccacgcgggaaacggtctgataagagacaccggcataactctgcgaca
tcgtataacgttactggtttcacattcaccaccctgaattgactctcttccgggcgctatcatg
ccataaccgcgaaagggttttgcgccattcgatggtgtccgggatctcgacgctctcccttatgcg
actcctgcattaggaagcagcccagtagtaggttgaggccgttgagcaccgcccgcgcaaggaa
tggtgcatgcaaggagatggcgcccaacagtcccccggccacggggcctgccaccataccacg
ccgaaacaagcgtcatgagcccgaagtggcgagcccgatcttccccatcggtgatgtcggcga
tataggcgccagcaaccgcacctgtggcgccggtgatgccggccacgatgcgtccggcgtagag
gatcgagatctcgatcccgcgaaattaatacgaactcactataggggaattgtgagcggataaca
attcccctctagaaataattttgtttaactttaagaaggagatataccatggagcataattccc
gtcgttccgcaaactatcagccaaacctgtggaatttcgaattcctgcaatccctggagaacga
cctgaaagtggaaaagctggaggagaaagcgaccaaactggaggaaagagttcgctgcagatc
aaccgtgtagacaccagccgctgtccctgctggagctgatcgacgatgtgcagcgccctgggtc
tgacctacaaatttgaaaaagacatcattaagccctggaaaacatcgtactgctggacgaaaa

Figure 68C

Nucleotide sequence of pDu50-4

caaaaagaacaaatctgacctgcacgcaaccgctctgtctttccgtctgctgcgtcagcacggt
ttcgagggtttctcaggatgtttttgagcggttcaaggataaagaagggtggtttcagcggtgaac
tgaaagggtgacgtccaaggcctgctgagcctgtatgaagcgtcttacctgggtttcgagggtga
gaacctgctggaggaggcgctaccttttccatcacccacctgaagaacaacctgaaagaaggc
attaataccaagggtgcagaacaagtgagccacgacctggaactgccatatacaccagcgtctgc
accgtctggaggcacgttggttcctggataaatacgaaccgaaagaaccgcatacaccagctgct
gctggagctggcgaagctggattttaacatggtacagacctgcaccagaaagagctgcaagat
ctgtcccgctggtggaccgagatggcctggctagcaaaactggattttgtacgcgaccgcctga
tggaagttaatttctgggcactgggtatggcgccagacccgcagtttggtgaatgtcgcaaagc
tgttactaaaatgtttggtctggtgacgatcatcgatgacgtgtatgacgtttatggcactctg
gacgaactgcaactgttcaccgatgctgtagagcgtgggacgttaacgctattaacacctgc
cggactatatgaaactgtgtttcctggcactgtacaacacccgttaacgacacgtcctattctat
tctgaaagagaaaaggtcataacaacctgtcctatctgacgaaaagctggcgtgaactgtgcaa
gcctttctgcaagaggcgaaatggtccaacaacaaaattatcccggttttctccaagtacctgg
aaaacgccagcgtttcctcctccggtgtagcgtgctggcgccgtcttacttttccgtatgcc
gcagcaggaagacatctccgaccacgcgctgcgttccctgaccgaactccatggtctggtgcgt
tctagctgcgttatcttccgcctgtgcaacgatctggccacctctgcggcgagctggaacgtg
gcgagactaccaattctatcattagctacatgcagaaaacgatggtaccagcgaggaacaggc
ccgcgaagaactgcgtaaactgatcgacgccgaatggaaaaagatgaatcgtgaacgcgttagc
gactccacctgctgcctaaagcgttcattggaaatcgacgttaacatggcacgtgtttccact
gcacctaccagtatggcgatggtctgggtcgcccagactacgcgactgaaaaccgcataaaact
gctgctgattgacctttcccgattaaccagctgatgtatgtctaaactgcagggatccgaattc
gagctccgtcgacaagcttgcgccgcac

Figure 69




Figure 70




Figure 71




Figure 72A

Nucleotide sequence of p9795 (SEQ ID NO:154)

tagaaaaactcatcgagcatcaaatgaaactgcaatttattcatatcaggattatcaat
accatatttttgaaaaagccgtttctgtaatgaaggagaaaaactcaccgaggcagttcc
ataggatggcaagatcctggtatcgggtctgcgattccgactcgtccaacatcaatacaa
cctattaatttcccctcgtcaaaaataagggttatcaagtgagaaatcaccatgagtgc
gactgaatccggtgagaatggcaaaagtttatgcatttctttccagacttggtcaacag
gccagccattacgctcgtcatcaaatcactcgcatcaaccaaaccgttattcattcgt
gattgcgcctgagcgcgagcgaaatacgcgatcgtgttaaaaggacaattacaaacagg
aatcgagtgcaaccggcgaggaacactgccagcgcatacaaatattttcacctgaat
caggatattcttctaataacctggaacgctgttttccggggatcgagtggtgagtaac
catgcatacatcaggagtagcgataaaatgcttgatggtcggaagtggcataaattccgt
cagccagtttagtctgaccatctcatctgtaacatcattggcaacgctaccttgccat
gtttcagaaacaactctggcgcatcgggcttccatacaagcgatagattgtcgacct
gattgcccgcacattatcgcgagcccatttatacccatataaatcagcatccatgttga
atttaatcgcgccctcgacgtttcccggttgaatatggctcatattcttcccttttcaat
attattgaagcatttatcagggttattgtctcatgagcggatacatatttgaatgtatt
tagaaaaataacaaataggggtcagtggtacaaccaattaaccaattctgaacattat
cgcgagcccatttataacctgaatatggctcataaacaccccttggttgccctggcggcag
agcgcggtggtcccacctgaccccatgccgaactcagaagtgaacgcgcgtagcgccga
tggtagtggtgggactccccatgagagagtagggaactgccaggcatcaaataaaacga
aaggctcagtcgaaagactgggcctttcgccgggctaattaggggggtgtcgcccttct
gattgacgctgcagttagacatacatcagctggttaatcgggaaagggtcaatcagcag
cagtttgatgcggttttcagtcgcgtagtctgggcgaccagaccatcgccatactggt
aggtgcagtgggaaacagtgccatgttaactgcgatttccatgaacgctttaggcagc
aggggtggagtgcgtaacgcgttcacgattcatcttttccattcggcgtcgatcagttt
acgcagttcttgcgggacctgttccctcgtggtaccatcgttttcgtgcagtagctaa
tgatagaattggtagtcctcgccaggttccagctccgcgcgagaggtggccagatcgttg
cacaggcggaagataacgcagctagaacgcaccagaccatggaagtcggtcagggaacg
cagcgcggtggtcgagatgtcttccctgctgctggcatacggaaaagtgaacggcgcca
gcagcgctacaccggaggaggaaacgctggcggttttccaggtacttgagaaaagccggg
ataattttgttggttggaacctttcgccctcttgcaaaaaggctttgcacagttcacgcca
gcttttcgtcagataggacaggttggttatgaccttctctttcagaatagaataggacg
tgctggttaacggtgttgtagcagtgccaggaaacacagtttcatatagtcgggcagggtg
ttaatagcgttaacgtcccagcgctctacagcatcgggaacagttgcagttcgtccag
agtgccataaacgtcatacacgtcatcgatgatcgtcaccagaccaaacttttagttaa
cagctttgcgacattcaccaaactgcgggtctggcgccataccagtgcccagaaataa
acttccatcaggcggtcgcgtaaaaaatccagtttgctagccaggcccatctcggtcca
ccagcgggacagatcttgagctctttctgggtgcagggtctgtaccatgttaaaatcca
gcttcgcccagctccagcagcagctgggtgatgcggttctttcggttcgtatttatccagg
aaccaacgtgcctccagacggtgcagacgctgggtgatatggcagttccaggggcgtggct
cacttggtctgcaaccttggtattaatgccttctttcaggttggttcttcaggtgggtga
tggaaaaaggtagcgccctcctccagcaggttctcacccctcgaaacccaggtgaacgct
tcatacaggtcagcaggccttggaacgtcacctttcagttcacccgctgaaaccaccttc
tttatccttgaaacgctcaaaaacatcctgagaaacctcgaaaccgtgctgacgcagca

Figure 72B

Nucleotide sequence of p9795

gacggaaagacagagcgggtgcgtgcaggtcagatttggtctttttgttttcgtccagc
agtacgatgttttccagggctttaatgatgtctttttcaaatttgtaggtcagacccag
gcgtgcacatcgtcgatcagctccagcaggacagcggctgggtgtctacacgggtga
tcatgcagcgaacttcttctccagtttggtcgctttctctccagcttttccacttc
aggtcgttctccagggattgcaggaattcgaaattccacaggttggtgatagttgc
ggaacgacgggaattatgctcggtaatctgagtaaattgagaagaggtcgacacatgt
tcagcgacaaggcgacacaaaaatttattctaaatgcataataaatactgataacatct
tatagtttggtattatattttgtattatcggtgacatgtataattttgatataaaaaact
gattttccctttattatttttcgagatttattttcttaattctctttaaaaaactagaaa
tattgtatatacaaaaaatcataaataatagatgaatagtttaattataggtgttcac
aatcgaaaaagcaacgtatcttattttaaagtgcgttgcttttttctcatttataaggtt
aaataattctcatatatacaagcaaaagtgcagggcgcccttaaataattctgacaaatgct
ctttccctaaactccccccataaaaaaacccgcgaagcgggtttttacgttatttgcg
gattaacgattactcgttatcagaaccgcccagggggcccagacttaagactggcgctc
gttttacaacacagaaagagtttgtagaaacgcaaaaaggccatccgtcaggggcttc
tgcttagtttgatgcctggcagttccctactctcgcttccgcttctcgtcactgac
tcgtgcgtcgggtcgttcgggtgcggcgagcggatcagctcactcaaaggcggtaac
acggttatccacagaatcaggggataacgcaggaaagaacatgtgagcaaaaggccagc
aaaaggccaggaaccgtaaaaaggcgcgttgctggcgtttttccataggctccgcccc
cctgacgagcatcacaaaaatcgacgctcaagtcagaggtggcgaaaccgcagaggact
ataaagataaccaggcgtttccccctggaagctccctcgtgcgtctcctgttccgacct
tgccgcttaccggataacctgtccgctttctcccttcgggaagcgtggcgctttctcat
agctcacgctgtaggtatctcagttcgggtgtaggtcgttcgctccaagctgggctgtgt
gcacgaaccccccggttcagcccagcgtgcgcttatccggttaactatcgtcttgagt
ccaaccggtaagacacgacttatcgccactggcagcagccactggtaacaggattagc
agagcgaggtatgtaggcgggtgctacagagttcttgaagtgggtgggctaactacggcta
cactagaagaacagtatattggtatctgcgctctgctgaagccagttaccttcggaaaaa
gagttggtagctcttgatccggcaaaaccacccgctggtagcgggtgggttttttgggt
tgcaagcagcagattacgcgcagaaaaaaaggatctcaagaagatcctttgatcttttc
tacgggggtctgacgctcagtggaacgacgcgcgtaactcacgttaagggttttgggt
catgagcttgcgccgtcccgtcaagtcagcgtaatgctctgcttt

Figure 73




Figure 74A

Nucleotide sequence of pTrcKudzu (SEQ ID NO:155)

gtttgacagcttatcatcgactgcacggtgcaccaatgcttctggcgtcaggcagccat
cggaagctgtggtatggctgtgcaggtcgtaaatcactgcataattcgtgtcgctcaag
gcgcactcccgttctggataatgtttttgcccgcacatcataacggttctggcaaata
ttctgaaatgagctgttgacaattaatcatccggtcgtataatgtgtggaattgtgag
cggataacaatttcacacaggaaacagcgccgctgagaaaaagcgaagcggcactgctc
tttaacaatttatcagacaatctgtgtgggcactcgaccggaattatcgattaacttta
ttattaaaaattaaagaggtatatattaatgtatcgattaaataaggaggaataaacca
tgtgtgcgacctcttctcaatttactcagattaccgagcataattcccgtcgttccgca
aactatcagccaaacctgtggaatttcgaattcctgcaatccctggagaacgacctgaa
agtggaaaagctggaggagaaagcgaccaaactggaggaagaagtccgtgcatgatca
accgtgtagacacccagccgctgtccctgctggagctgatcgacgatgtgcagcgctg
ggtctgacctacaaatttgaaaaagacatcattaaagccctggaaaacatcgactgct
ggacgaaaacaaaaagaacaaatctgacctgcacgcaaccgctctgtcttccgtctgc
tgcgtcagcacggtttcgagggtttctcaggatgtttttgagcgtttcaaggataaagaa
ggtggtttcagcgggtgaactgaaaggtgacgtccaaggcctgctgagcctgtatgaagc
gtcttacctgggtttcgagggtgagaacctgctggaggaggcgcgtaccttttccatca
cccacctgaagaacaacctgaaagaaggcattaataccaaggttgagaacaagtgagc
cacgccctggaactgccatatcaccagcgtctgcaccgtctggaggcacgttggttct
ggataaatacgaaccgaaagaaccgcatcaccagctgctgctggagctggcgaagctgg
attttaacatggtacagaccctgcaccagaaagagctgcaagatctgtcccgtggtgg
accgagatgggcctggctagcaaaactggattttgtacgcgaccgctgatggaagtta
tttctgggcaactgggtatggcgccagaccgcagtttggtgaatgtcgaaagctgtta
ctaaaatgtttggtctggtgacgatcatcgatgacgtgtatgacgtttatggcactctg
gacgaaactgcaactgttcaccgatgctgtagagcgtgggacgttaacgctattaacac
cctgccggactatatgaaactgtgtttcctggcactgtacaacaccgttaacgacacgt
cctattctattctgaaagagaaaggtcataacaacctgtcctatctgacgaaaagctgg
cgtgaaactgtgcaaagcctttctgcaagaggcgaaatggtccaacaacaaaattatccc
ggctttctccaagtacctggaaaacgccagcgtttcctcctccggtgtagcgtgctgg
cgccgtcttacttttccgtatgccagcagcaggaagacatctccgaccacgcgtgcgt
tccctgaccgacttccatggtctggtgcttctagctgcgttatcttccgcctgtgcaa
cgatctggccacctctgcggcgagctggaacgtggcgagactaccaattctatcatta
gctacatgcacgaaaacgatggtaccagcgaggaacaggcccgcgaagaactgcgtaaa
ctgatcgacgcccgaatggaaaaagatgaatcgtgaacgcgttagcgactccacctgct
gcctaaagcgttcatggaaatcgagttaacatggcacgtgttcccactgcacctacc
agtatggcgatggtctgggtcgcccagactacgcgactgaaaaccgcatcaaactgctg
ctgattgaccctttcccgattaacagctgatgtatgtctaactgcagctggtaccata
tgggaattcgaagctttctagaacaaaaactcatctcagaagaggatctgaatagcgcc
gtcgaccatcatcatcatcattgagtttaaaccggtctccagcttggtgttttggc
ggatgagagaagattttcagcctgatacagattaaatcagaacgcagaagcggctgtat
aaaacagaatttgccctggcggcagtagcgcggtggtcccacctgaccccatgccgaact
cagaagtgaacgccgtagcgccgatggtagtgtgggtctccccatgcgagagtaggg
aactgccaggcatcaaataaaacgaaaggctcagtcgaaagactgggcctttcgtttta
tctgttgtttgtcggtgaacgctctcctgagtaggacaaatccgcccggagcggatttg

Figure 74B

Nucleotide sequence of pTrcKudzu

aacgttgcgaagcaacggcccgagggtggcgggcaggacgcccgcataaaactgccag
gcatcaaattaagcagaaggccatcctgacggatggcctttttgcgtttctacaaactc
tttttggtttatttttctaaatacattcaaataatgtatccgctcatgagacaataaccct
gataaatgcttcaataaatattgaaaaaggaagagtatgagtattcaacatttccgtgtc
gcccttattcccttttttgcggcattttgccttccgtgttttgcctcaccagaaaacgct
ggtgaaagtaaaagatgctgaagatcagttgggtgcacgagtgggttacatcgaactgg
atctcaacagcggtaagatccttgagagttttcgccccgaagaacgtttccaatgatg
agcacttttaagttctgctatgtggcgcggtattatcccggtgttgacgcgggcaaga
gcaactcggtcgcccgcatacactattctcagaatgacttgggtgagtactcaccagtca
cagaaaagcatcttacggatggcatgacagtaagagaattatgcagtgtgccataacc
atgagtataacactgcggccaacttacttctgacaacgatcggaggaccgaaggagct
aacgcttttttgacacaacatgggggatcatgtaactcgcttgatcgttgggaaccgg
agctgaatgaagccataccaaacgacgagcgtgacaccacgatgcctgtagcaatggca
acaacgttgcgcaaactattaactggcgaactacttactctagcttcccggaacaatt
aatagactggatggaggcggataaagttgcaggaccacttctgcgctcgcccttccgg
ctggctggtttattgctgataaatctggagccggtgagcgtgggtctcgcggtatcatt
gcagcaactggggccagatggttaagccctcccgatcgtagtattctacacgacggggag
tcaggcaactatggatgaacgaaatagacagatcgctgagatagggtgectcactgatta
agcattggtaactgtcagaccaagtttactcatatatacttttagattgatttaaaactt
catttttaatttaaaaggatctaggtgaagatcctttttgataatctcatgacccaaat
cccttaacgtgagttttcgttccactgagcgtcagaccccgtagaaaagatcaaaggat
cttcttgagatccttttttctgcgcgtaaatctgctgcttgcaaacaaaaaaaccaccg
ctaccagcgggtggtttgtttgcgggatcaagagctaccaactcttttccgaaggtaac
tggttccagcagagcgcagataccaaataactgtccttctagtgtagcctagttaggcc
accaactcaagaactctgtagcaccgcctacatacctcgctctgctaactcctgttacca
gtggctgctgccagtggcgataagtcgtgtcttaccgggttggaactcaagacgatagtt
accggataaaggcgcagcggctcggtgtgaacggggggttcgtgcacacagcccagcttg
agcgaacgacctacaccgaactgagatacctacagcgtgagctatgagaaagcggcacg
cttcccgaaggagagaaaggcggacaggtatccggttaagcggcagggtcggaacaggaga
gcgcaogaggggagcttccagggggaacgcctggtatctttatagtcctgtcgggttc
gccacctctgacttgagcgtcgatttttgtgatgctcgtcaggggggaggagcctatgg
aaaaacgccagcaacgcggcctttttacggttcctggccttttgctggccttttgctca
catgttctttcctgcgttatccctgattctgtggataaacgattaccgcctttgagt
gagctgataccgctcgccgcagccgaacgaccgagcgcagcagtcagtgagcaggaa
gcggaagagcgcctgatgcggtattttctccttacgcatctgtgcggtatttcacaccg
catatggtgcactctcagtacaatctgctctgatgccgcatagttaagccagtatacac
tccgctatcgctacgtgactgggtcatggctgcgccccgacaccgcgaacaccgcgtg
acgcgcctgacgggcttgtctgctcccggcatccgcttacagacaagctgtgaccgtc
tccgggagctgcatgtgtcagaggttttcaccgtcatcaccgaaacgcgcgaggcagca
gatcaattcgcgcggaaggcgaagcggcatgcatttacgttgacaccatcgaatgggtg
caaaacctttcgcggtatggcatgatagcggccggaagagagtcaattcaggggtggtga
atgtgaaaccagtaacgttatacgatgtcgagagtatgccggtgtctcttatcagacc
gtttcccgctggtgaaccaggccagccacgttctgcgaaaacgcgggaaaaagtga

Figure 74C

Nucleotide sequence of pTrcKudzu

agcggcgatggcggagctgaattacattcccaaccgcgtggcacaacaactggcgggca
aacagtcgttgctgattggcgttgccacctccagctctggccctgcacgcgccgtcgcaa
attgtcgcggcgattaaatctcgcgccgatcaactgggtgccagcgtggtggtgtcgat
ggtagaacgaagcggcgtcgaagcctgtaaagcggcgggtgcacaatcttctcgcgcaac
gcgtcagtgggctgatcattaactatccgctggatgaccaggatgccattgctgtggaa
gctgcctgcaactaatgttccggcgttatttcttgatgtctctgaccagacacccatcaa
cagtattatttctcccatgaagacggtagcgcgactgggcgtggagcatctggtcgcat
tgggtcaccagcaaatcgcgctgttagcgggcccattaaagtctgtctcggcgcgtctg
cgtctggctggctggcataaatatctcactcgcaatcaaattcagccgatagcggaacg
ggaaggcgactggagtgccatgtccggttttcaacaaaccatgcaaatgctgaatgagg
gcatcgttcccactgcgatgctggttgccaacgatcagatggcgctgggcgcaatgcgc
gccattaccgagtcgggctgcgcgttggtgcggaatatctcggtagtgggatacgacga
taccgaagacagctcatgttatatccgcgcgtcaaccaccatcaaacaggattttcgcc
tgctggggcaaaccagcgtggaccgcttgctgcaactctctcagggccaggcgggtgaag
ggcaatcagctgttgcccgctctcactggtgaaaagaaaaaccaccctggcgcccaatac
gcaaaccgcctctccccgcgcgttggcggattcattaatgcagctggcacgacaggttt
cccgactggaaagcgggcagtgagcgcaacgcaattaatgtgagttagcgcgaattgat
ctg

Figure 75




Figure 76A

Nucleotide sequence of pMAL-C4X (SEQ ID NO:156)

ccgacaccatcgaatggtgcaaaacctttcgcggtatggcatgatagcgccccgaagag
agtcaattcaggggtggtgaatgtgaaaccagtaacgttatacgatgtcgagagatgc
cggtgtctcttatcagaccgtttcccgcgtggtgaaccaggccagccacgtttctgcga
aaacgcgggaaaaagtggaagcggcgatggcggagctgaattacattcccaaccgcgtg
gcacaacaactggcgggcaaacagtcgttgctgattggcgttgccacctccagctctggc
cctgcacgcgcgtcgcaaattgtcgcggcgattaaatctcgcgccgatcaactgggtg
ccagcgtggtggtgctgatggtagaacgaagcggcgctgaagcctgtaaagcggcggtg
cacaatcttctcgcgcaacgcgtcagtggtgatcattaaactatccgctggtgacca
ggatgccattgctgtggaagctgctgcaactaatgttccggcggttatttcttgatgtct
ctgaccagacacccatcaacagttatttttctccatgaagacggtagcgactgggc
gtggagcatctggtcgattgggtcaccagcaaactcgcgctgttagcgggcccattaa
ttctgtctcggcgcgtctgcgtctggtggctggcataaatatctcactcgcaatcaaa
ttcagccgatagcggaaacgggaaggcgaactggagtgccatgtccggttttcaacaaacc
atgcaaatgctgaatgagggcatcggttcccactgcgatgctggttgccaacgatcagat
ggcgctgggcgcaatgcgcgccattaccgagtcggggctgcgcgttggtgcggatatct
cggtagtgggatacgacgataccgaagacagctcatgttatatcccgcgttaaccacc
atcaaacaggattttcgccctgctggggcaaacaccagcgtggaccgcttgctgcaactctc
tcaggggccaggcgggtgaagggcaatcagctgttgcccgctctcactggtgaaaagaaaa
ccacctggcgcccaatacgcaaacgcctctcccgcgcgttgcccgattcattaatg
cagctggcacgacaggtttcccgaactggaaagcgggcagtgagcgcaacgcaattaatg
taagttagctcactcattaggcacaattctcatgtttgacagcttatcatcgactgcac
ggtgcaccaatgcttctggcgtcaggcagccatcggaagctgtggtatggctgtgcagg
tcgtaaatcactgcataattcgtgtcgctcaaggcgcactccggttctggataatgttt
tttgcgcgcgacatcataacggttctggcaaatattctgaaatgagctgttgacaattaa
tcacggtcgtataatgtgtggaattgtgagcggataacaatttcacacaggaaacag
ccagtcggttttaggtgttttcacgagcacttcaccaacaaggaccatagcatatgaaaa
tcgaagaaggtaaactggtaatctggattaacggcgataaaggctataacggtctcgct
gaagtgcgtgaagaaattcgagaaagataccggaattaaagtaccggttgagcatccgga
taaactggaagagaaattcccacaggttgccggcaactggcgatggccctgacattatct
tctggggcacacgacgcgtttggtggctacgctcaatctggcctgttggtgaaatcacc
ccggacaaagcgttccaggacaagctgtatccgtttacctgggatgccgtacgttacaa
cggcaagctgatttgcttaccgatcgctgttgaaagcgttatcgctgatttataacaaag
atctgctgccgaacccgcacaaaacctgggaagagatcccggcgctggataaagaactg
aaagcgaaaggtgaagagcgcgtgatgttcaacctgcaagaaccgtacttcacctggcc
gctgattgctgctgacgggggttatgcgttcaagtatgaaaacggcaagtacgacatta
aagacgtgggcgtggataacgctggcgcgaaagcgggtctgacctcctggttgacctg
attaaaaacaaacacatgaatgcagacaccgattactccatcgagaagctgcctttaa
taaaggcgaaacagcgatgacctcaacggcccgtgggcatggtccaacatcgacacca
gcaaagtgaattatggtgtaacggtactgccgaccttcaagggtcaaccatccaaaccg
ttcggtggcgtgctgagcgcaggtattaacgcgcgacgtccgaacaaagagctggcaaa
agagttcctcgaaaactatctgctgactgatgaaggtctggaagcgttaataaagaca
aacgcgtgggtgccgtagcgtgaagcttacgaggaagagttggtgaaagatccgcgg
attgccgccactatggaaaacgcccagaaaggtgaaatcatgccgaacatcccgcagat

Figure 76B

Nucleotide sequence of pMAL-C4X

gtccgctttctggtatgccgtgcgtactgcggtgatcaacgccgccagcggcgtcgcagac
ctgtcgatgaagccctgaaagacgcgcagactaattcgagctcgaacaacaacaacaat
aacaataacaacaacctcgggatcgaggggaaggatttcagaattcggatcctctagagt
cgacctgcaggcaagcttggcactggccgtcgttttacaacgtcgtgactgggaaaacc
ctggcgttacccaacttaatcgcttgcagcacatcccccttcgccagctggcgtaat
agcgaagaggccgcaccgatcgcccttcccaacagttgcgcagcctgaatggcgaatg
gcagcttggctgttttggcggatgagataagattttcagcctgatacagattaaatcag
aacgcgaagcggctctgataaaacagaatttgcctggcggcagtagcgcggtggtccca
cctgaccccatgccgaactcagaagtgaacgccgtagcgccgatggtagtgtggggtc
tccccatgcgagagtagggaactgccaggcatcaaataaaacgaaaggctcagtcgaaa
gactgggcctttcgttttatctgttgttgcggtgaacgctctcctgagtaggacaaa
tccgcccgggagcggatttgaacgttgcgaagcaacggcccgagggtggcgggcaggac
gcccgcataaactgccaggcatcaaattaagcagaaggccatcctgacggatggcctt
tttgcgtttctacaaactcttttgtttatttttctaaatacattcaaataatgtatccg
ctcatgagacaataaacctgataaatgcttcaataaatattgaaaaaggaagagtatgag
tattcaacatttcggtgtcgcccttattccctttttgcggcattttgccttcctgttt
ttgctcaccagaaaacgctggtgaaagtaaaagatgctgaagatcagttgggtgcacga
gtgggttacatcgaactggatctcaacagcggtaagatccttgagagttttcgccccga
agaacgttctccaatgatgagcacttttaaagttctgctatgtggcgcggtattatccc
gtgttgacgcgggcaagagcaactcggtcgcgcatacactattctcagaatgacttg
gttgagtactcaccagtacagaaaagcatcttacggatggcatgacagtaagagaatt
atgcagtgctgccataaacatgagtataacactgcggccaacttacttctgacaacga
tcggaggaccgaaggagctaacgcgttttttgcacaacatgggggatcatgtaactcgc
cttgatcgttgggaaccggagctgaatgaagccataccaaacgacgagcgtgacaccac
gatgcctgtagcaatggcaacaacgttgcgcaaacatttaactggcgaactacttactc
tagcttcccggcaacaattaatagactggatggaggcggataaagttgcaggaccactt
ctgcgctcggcccttccggtggtgtttattgctgataaatctggagccggtgagcg
tgggtctcgcggtatcattgcagcactggggccagatggtgaagccctccgctatcgtag
ttatctacacgacggggagtcaggcaactatggatgaacgaaatagacagatcgctgag
ataggtgcctcactgattaagcattggttaactgtcagaccaagtttactcatatatact
ttagattgatttaccocgggttgataatcagaaaagccccaaaaacaggaagattgtata
agcaaatatttaaatgttaaacttaataattttgttaaaattcgcgttaaattttgtt
aaatcagctcattttttaaccaataggccgaaatcggaacaaatcccttataaatcaaaa
gaatagaccgagatagggttgagtggttccagtttggacaagagtcactattaaa
gaacgtggactccaacgtcaaagggcgaaaaaccgtctatcagggcgatggcccactac
gtgaaccatcacccaaatcaagtttttggggcgaggtgccgtaaagcactaaatcgg
aaccttaaagggagcccccgatttagagcttgacggggaaagccggcgaacgtggcgag
aaaggaaggggaagaaagcgaaaggagcgggcgctagggcgctggcaagtgtagcggca
cgctgcgcgtaaccaccacaccgcgcgcttaatgcgcgcgtacagggcgcgtaaaag
gatctaggtgaagatcctttttgataatctcatgacaaaatcccttaacgtgagtttt
cgttccactgagcgtcagaccccgtagaaaagatcaaaggatcttcttgagatcctttt
tttctgcgcgtaatctgctgcttgcaacaaaaaaaccaccgctaccagcgggtggttg
tttgccggatcaagagctaccaactctttttccgaaggttaactggcttcagcagagcgc

Figure 76C

Nucleotide sequence of pMAL-C4X

agataccaaatactgtccttctagtgtagccgtagtttaggccaccacttcaagaactct
gtagcaccgcctacatacctcgctctgctaatacctgttaccagtggctgctgccagtgg
cgataagtctgtcttaccgggttggaactcaagacgatagttaccggataaaggcgcagc
ggtcgggctgaacgggggggttcgtgcacacagcccagcttggagcgaacgacctacacc
gaactgagatacctacagcgtgagctatgagaaagcgccacgcttcccgaaggagaaa
ggcggacaggtatccggt aagcggcaggggtcggaacaggagagcgcacgaggagcttc
cagggggaaacgcctggatctttatagtccctgtcgggtttcgccacctctgacttgag
cgtcgatttttgtgatgctcgtcaggggggaggagcctatggaaaaacgccagcaacgc
ggcctttttacgggttcctggccttttgcctggccttttgcctcacatgttcttctcgt
tatcccctgattctgtggataaacgctattaccgcctttgagtgagctgataccgctcgc
cgcagccgaacgaccgagcgcagcagctcagtgagcaggaagcgggaagagcgctgat
gcggtattttctccttacgcctctgtgcgggtatttcacaccgcataatggtgcactct
cagtacaatctgctctgatgccgcatagttaagccagtatataactccgctatcgctacg
tgactgggtcatggtgcgccccgacacccgccaaacacccgctgacgcgcctgacggg
cttgtctgctcccggcatccgcttacagacaagctgtgaccgtctccgggagctgcatg
tgtcagagggtttcacccgtcatcaccgaaacgcgcgaggcagctgcggtaagctcatc
agcgtggctcgtgcagcgtatcacagatgtctgectgttcatccgcgtccagctcgttga
gtttctccagaagcgttaatgtctggttcttgataaagcggggccatgttaagggcgggt
ttttcctgtttgggtcactgatgcctccgtgtaaggggatttctgttcatgggggtaat
gataccgatgaaacgagagaggatgctcacgatacgggttactgatgatgaacatgcc
ggttactggaacgttgtgagggtaacaactggcggatggatgcggcgggaccagaga
aaaatactcagggtaatgccagcgttcgttaatacagatgtagggtgttccacaggg
tagccagcagcatcctgcgatgcagatccggaacataatggtgcagggcgctgacttcc
gcgtttccagactttacgaaacacggaaacgaagaccattcatgttgttgcctcaggtc
gcagacgttttgcagcagcagctcgttcacgttcgctcgcgtatcggtgattcattctg
ctaaccagtaaggcaacccccgccagcctagccgggtcctcaacgacaggagcacgatca
tgcgcacccgtggccaggacccaacgctgcccgaaatt

Figure 77




Figure 78A

Nucleotide sequence of pMAL-C4X Kudzu (SEQ ID NO:157)

cgacacaccatcgaaatggtgcaaaaccttttcgcggtatggcatgatagcgccgcggaagag
agtcaattcaggggtggtgaatgtgaaaccagtaacgttatacgatgtcgagagtatgc
cgggtgtctcttatcagaccgtttcccgcggtggtgaaccaggccagccacgtttctgcga
aaacgcgggaaaaaagtggaaagcggcgatggcggagctgaattacattcccaaccgcgtg
gcacaacaactggcgggcaaacagtcgttgcgtattggcgttggcacctccagtcggc
cctgcacgcgcgtcgcaaattgtcgcgggcattaaatctcgcgccgatcaactgggtg
ccagcgtggtggtgtcgatggtagaacgaagcggcgtcgaagcctgtaaagcggcggtg
cacaatctctcgcgcaacgcgtcagtgggctgatcattaactatccgctggatgacc
ggatgccattgctgtggaagctgctgcactaatgttcggcgcttatttcttgatgtc
ctgaccagacacccatcaacagtatatttttctcccatgaagacggtagcgactgggc
gtggagcatctggtgcattgggtcaccagcaaatcgcgctgttagcgggcccattaaag
ttctgtctcggcgcgctctcgctctggctggctggcataaatatctcactcgcaatcaaa
ttcagccgatagcggaaacgggaaggcgactggagtcccatgtcgggttttcaacaaacc
atgcaaatgctgaatgagggcatcgttcccactgcgatgctggttggcaacgatcagat
ggcgctgggcgcaatgcgcgccattaccagctccgggctgcgcgttggtgcggatatct
cggtagtgggatacgacgataccgaagacagctcatgttatatccgcgcgttaaccacc
atcaaacaggatttttcgcctgctggggcaaacacgcgtggaccgcttgctgcaactctc
tcagggccaggcgggtgaagggcaatcagctgttgcccgtctcactggtgaaaagaaaa
ccacctggcgcccaatacgcaaaccgcctctccccgcgcgttggccgattcattaatg
cagctggcacgacaggtttcccgactggaaagcgggcagtgagcgcaacgcaattaatg
taagttagctcactcattaggcacaatttctcatgtttgacagcttatcatcgactgcac
ggtgcaccaatgcttctggcgtcaggcagccatcggaagctgtggtatggctgtgcagg
tcgtaaatcactgcataattcgtgtcgctcaaggcgcaactccggttctggataatgttt
tttgcgcgcacatcataacggttctggcaaatattctgaaatgagctgttgacaattaa
tcatcggctcgtataatgtgtggaattgtgagcggataacaatttcacacaggaacag
ccagtcggtttaggtgttttcagcagcacttcaccaacaaggaccatagcatatgaaa
tcgaagaaggtaaactggtaatctggattaacggcgataaaggctataacggctctcgct
gaagtccgttagaaattcgagaaagataccggaattaaagtaccggttagcatccgga
taaactggaagagaaattcccacaggttgcggcaactggcgatggccctgacattatct
tctgggcacacgaccgctttggtggctacgctcaatctggcctgttggctgaaatcacc
ccggacaaagcgttcaggacaagctgtatccgtttacctgggatgccgtacgttacaa
cggcaagctgattgcttaccgatcgctgttgaagcgttatcgctgatttataacaaag
atctgctgccgaaccgcgcaaaaacctgggaagagatccgcgcgtggataaagaactg
aaagcgaaaggttagagcgcgtgatgttcaacctgcaagaaccgtacttcacctggc
gctgattgctgctgacgggggttatgcgttcaagtatgaaaacggcaagtagcatata
aagacgtggcggtggataacgctggcgcgaaagcgggtctgaccttccctggttagcctg
attaaaaacaaacacatgaatgcagacaccgattactccatcgagaagctgcctttaa
taaggcgaaacagcgatgaccatcaacggcccgtgggcatggtccaacatcgacacca
gcaaagtgaattatggtgtaacggtagtgcggaccttcaagggtcaaccatccaacccg
ttcgttggcgtgctgagcgcaggtattaacgcgccagtcgaacaaagagctggcaaa
agagttcctcgaaaactatctgctgactgatgaaggctggaagcgggttaataaagaca
aacccgctgggtgccgtagcgcgtgaagtcttacgaggaagagttggtgaaagatccgcgg
attgccgccactatggaaaacgccagaaaggtgaaatcatgccgaacatccgcgat

Figure 78B

Nucleotide sequence of pMAL-C4X Kudzu

gtccgctttctggtatgccgtgcgtactgcggtgatcaacgccgccagcggtcgtcaga
ctgtcgatgaagccctgaaagacgcgcagactaattcgagctcgaacaacaacaacat
aacaataacaacaacctcgggatcgagggaaggatttcagaattctgtgcgacctcttc
tcaatttactcagattaccgagcataattcccgtcgttcgcgaaactatcagccaaacc
tgtggaatttcgaattcctgcaatccctggagaacgacctgaaagtggaaaagctggag
gagaaagcgaccaaaactggaggaagaagttcgtgcatgatcaaccgtgtagacacca
gccgctgtccctgctggagctgatcgacgatgtgcagcgccctgggtctgacctacaaat
ttgaaaaagacatcattaaagccctggaaaacatcgtactgctggacgaaaacaaaaag
aacaatctgacctgcacgcaaccgctctgtctttccgtctgctgcgtcagcacggttt
cgaggtttctcaggatgtttttgagcggtttcaaggataaagaagggtggtttcagcggtg
aactgaaaggtgacgtccaaggcctgctgagcctgtatgaagcgtcttacctgggtttc
gaggggtgagaacctgctggaggaggcgctaccttttccatcacccacctgaagaacaa
cctgaaagaaggcattaataccäagggtgcagaacaagtgagccacgccctggaactgc
catatcaccagcgtctgcaccgtctggaggcacgttggttcctggataaatacgaaccg
aaagaaccgcatcaccagctgctgctggagctggcgaagctggattttaacatggtaca
gaccctgcaccagaaagagctgcaagatctgtcccgtcgttgaccgagatgggcctgg
ctagcaaaactggattttgtacgcgaccgcctgatggaagtttatttctgggcaactgggt
atggcgccagaccgcagtttggtgaatgtcgaaagctgttactaaaatgtttggtct
ggtgacgatcatcgatgacgtgtatgacgtttatggcactctggacgaactgcaactgt
tcaccgatgctgtagagcgctgggacgttaacgctattaacaccctgccggactatatg
aaactgtgtttcctggcactgtacaacaccgttaacgacacgtcctattctattctgaa
agagaaagggtcataacaacctgtcctatctgacgaaaagctggcgtgaactgtgcaaag
cctttctgcaagaggcgaaatggtccaacaacaaaattatcccggctttctccaagtac
ctggaaaacgccagcgtttcctcctccggtgtagcgtgctggcgccgtcttacttttc
cgtatgccagcagcaggaagacatctccgaccacgcgctgcgttcctgaccgacttcc
atggtctggtgcgttctagctgcgttatcttccgcctgtgcaacgatctggccacctct
gcgccggagctggaacgtggcgagactaccaattctatcattagctacatgcacgaaaa
cgatggtaccagcgaggaacaggcccgcaagaactgcgtaaactgatcgacgccgaat
ggaaaaagatgaatcgtgaacgcgttagcgactccaccctgctgcctaaagcgttcatg
gaaatcgaggttaacatggcacgtgtttccactgcacctaccagtatggcgatggtct
gggtcgcccagactacgcgactgaaaaccgcatcaaactgctgctgattgacctttcc
cgattaaaccagctgatgtatgtctaagcttggcactggccgtcgttttacaacgtcgtg
actgggaaaaccctggcggttaccctaatcgccttgcagcacatccccctttcgcc
agctggcgtaatagcgaagaggcccgaccgatcgcccttcccaacagttgcgcagcct
gaatggcgaatggcagcttggtgttttggcggtgagataagattttcagcctgatac
agattaaatcagaacgcagaagcgggtctgataaaacagaatttgccctggcggcagtagc
gcggtggtcccacctgaccccatgccgaactcagaagtgaacgccgtagcgccgatgg
tagtgtgggtctccccatgcgagagtagggaactgccaggcatcaaataaaacgaaag
gctcagtcgaaagactgggcctttcggtttatctgttgtttgtcggtgaacgctctcct
gagtaggacaaatccgcccgggagcggatttgaaacgttgcaagcaacggcccgagggt
ggcggggcaggacgcccgcataaactgccaggcatcaaattaagcagaaggccatcctg
acggatggcctttttgcgtttctacaaactctttttgtttatttttctaaatacattca
aatatgtatccgctcatgagacaataaccctgataaatgcttcaataatattgaaaaag

Figure 78C

Nucleotide sequence of pMAL-C4X Kudzu

gaagagtatgagtattcaacatttccgtgtcgcccttattcccttttttgcggcatttt
gccttcctgttttgcctcaccagaaaacgctggtgaaagtaaaagatgctgaagatcag
ttgggtgcacgagtggttacatcgaactggatctcaacagcggtaagatccttgagag
ttttcgccccgaagaacgttctccaatgatgagcacttttaaagtctgctatgtggcg
cggattatcccgtgttgacgcccgggcaagagcaactcggcgccgcatacactattct
cagaatgacttggttgagtactcaccagtacagaaaagcatcttacggatggcatgac
agtaagagaattatgcagtgtgcccataaccatgagtataaactgcggccaacttac
ttctgacaacgatcggaggaccgaaggagctaaccgcttttttgacaacatgggggat
catgtaactcgcccttgatcgttgggaaccggagctgaatgaagccataccaaacgacga
gcgtgacaccacgatgcctgtagcaatggcaacaacgttgcgcaaactattaactggcg
aactacttactctagcttcccggcaacaattaatagactggatggaggcggataaagtt
gcaggaccacttctgcgtcgcccttccggctggctggtttattgctgataaatctgg
agccggtgagcgtgggtctcgcggtatcattgcagcactggggccagatggtgaagccct
cccgtatcgtagttatctacacgacggggagtcaggcaactatggatgaacgaaataga
cagatcgctgagataggtgcctcactgattaagcattggttaactgtcagaccaagttta
ctcatatatacttttagattgatttaccgcggttgataatcagaaaagccccaaaaacag
gaagattgtataagcaaatatttaaattgtaaacgttaatattttgttaaaattcgct
taaatttttgttaaatcagctcattttttaaccaataggccgaaatcggcaaaatccct
tataaatcaaaagaatagaccgagataggggttgagtgtgttccagtttggacaagag
tccactattaaagaacgtggactccaacgtcaaagggcgaaaaaccgtctatcaggggcg
atggccactacgtgaaccatcacccaaatcaagtttttggggcgaggtgccgtaaa
gcactaaatcggaaccctaaagggagcccccgatttagagcttgacggggaaagccggc
gaacgtggcgagaaaggaagggaagaaagcgaaaggagcgggcgctagggcgctggcaa
gtgtagcgggtcacgctgcgcgttaaccaccacacccgcgcgcttaatgcgcgcgtacag
ggcgcgtaaaaggatctaggtgaagatccttttgataatctcatgacaaaaatccctt
aacgtgagttttcggtccactgagcgtcagaccccgtagaaaagatcaaaggatcttct
tgagatccttttttctgcgcgtaatctgctgcttgcaaacaaaaaaaccaccgctacc
agcgggtggtttgtttgocggatcaagagctaccaactcttttccgaaggtaactggct
tcagcagagcgcagataccaaatactgtccttctagttagccgtagttaggccaccac
ttcaagaactctgtagcacccgctacatacctcgctctgctaatcctgttaccagtggc
tgctgccagtggcgataagtcgtgtcttaccgggttggaactcaagacgatagttaccgg
ataaggcgagcgggtcgggctgaacgggggggttcgtgcacacagcccagcttggaagcga
acgacctacaccgaactgagatacctacagcgtgagctatgagaaagcgccacgcttcc
cgaagggaagaaaggcggacaggtatccggtgaagcggcagggtcggaacaggagagcgca
cgaggagccttccagggggaaacgcctggtatctttatagtcctgtcgggtttcgccac
ctctgacttgagcgtcgatttttgtgatgctcgtcagggggggcgagcctatggaaaaa
cgccagcaacgcggcctttttacgggttctggccttttgctggccttttgctcacatgt
tctttcctgcgttatcccctgattctgtggataaccgattaccgcctttgagtgaact
gataccgctcgccgcagccgaacgaccgagcgcagcagtcagtgaagcaggaagcgga
agagcgctgatgcggatattttctccttacgcctctgtgcggatattcacaccgcatat
atggtgcaactctcagtacaatctgctctgatgccgcatagttgaagccagtatacactcc
gctatcgctacgtgactgggtcatggctgcgccccgacacccgccaacacccgctgacg
cgccctgacgggcttgtctgctcccggcatccgcttacagacaagctgtgaccgtctcc

Figure 78D

Nucleotide sequence of pMAL-C4X Kudzu

gggagctgcatgtgtcagagggtttcaccgatcatcaccgaaacgcgcgaggcagctgcg
gtaaagctcatcagcgtggtcgtgcagcgattcacagatgtctgcctgttcacccgct
ccagctcgttgagtttctccagaagcgtaaatgtctggcttctgataaagcgggccatg
ttaaggcggttttttccctgtttggtcactgatgcctccgtgtaagggggatttctgtt
catgggggtaatgataccgatgaaacgagagaggatgctcacgatacgggttactgatg
atgaacatgcccggttactggaacgttgtgagggtaaacaactggcggatggatgcgg
cgggaccagagaaaaatcactcagggccaatgccagcgttcgttaatacagatgtagg
tgttccacagggtagccagcagcatcctgcgatgcagatccggaacataatggtgcagg
gcgctgacttccgcgtttccagactttacgaaacacggaaccgaagaccattcatgtt
gttgctcaggtcgcagacgttttgcagcagcagtcgcttcacgttcgctcgcgtatcgg
tgattcattctgctaaccagtaaggcaaccccgcagcctagccgggtcctcaacgaca
ggagcacgatcatgcgcacccgtggccaggacccaacgctgcccgaaatt

Figure 79A




Figure 79B




Figure 80

Amino acid sequence of *P. tremuloides* IspS (SEQ ID NO:158)

MRCSVSTENVSFSETETETRRSANYEPNSWDYDYLLSSDTDESIEVHKDKAKKLEAEVR
REINNEKAEFLTLLELIDNVQRLGLGYRFESDIRRALDRFVSSGGFDGVTKTSLHGTA
SFRLLRQHGFVVSQEAESGFKDQNGNFLENLKEDIKAILSLYEASFLALEGENILDEAK
VFAISHLKELSEEKIGKELAEQVSHALELPLHRRRTQRLEAVWSIEAYRKEDANQVLE
LAILDYNNMIQSVYQRDLRETSRWRRVGLATKLHFARDRLIESFYWAVGVAFEPQYSDC
RNSVAKMFSFVTIIDDYDVYGTLDLELEFTDAVERWDVNAINDLDPDYMKLCFLALYNT
INEIAYDNLKDKGENILPYLTAWADLCNAFLQEAKWLYNKSTPTFDDYFGNAWKSSSG
PLQLIFAYFAVVQNIKKEEIENLQKYHDIISRPSHIFRLCNDLASASAEIARGETANSV
SCYMRKGI SEELATESVMNLIDETWKKMNKEKLGGS LFAKPFVETAINLARQSHCTYH
NGDAHTSPDELTRKRVLSVITEPILPFER

Figure 81A

Nucleotide sequence of *P. tremuloides* pET24a (SEQ ID NO:159)

tggcgaatgggacgcgcctgtagcggcgcattaaagcgcggcgggtgtgggtggttacgc
gcagcgtgacgcgtacacttgccagcgcctagcgcgcgcctcctttcgctttcttcct
tcctttctcgccacgttcgcggctttccccgtcaagctctaaatcgggggctcccttt
agggttccgatttagtgctttacggcacctcgaccccaaaaaacttgattagggatgag
gttcacgtagtgggccaatcgccctgatagacgggttttcgcctttgacgttgaggtcc
acgttctttaatagtggactcttggtccaaactggaacaacactcaaccctatctcgg
ctattcttttgattataaagggttttgccgatttcggcctattgggttaaaaaatgagc
tgatttaacaaaaatttaacgcgaattttaacaaaatattaacgtttacaatttcagg
ggcacttttcggggaatgtgcgcggaacccctatttggtttattttctaaatacatc
aaatatgtatccgctcatgaattaatcttagaaaaactcatcgagcatcaaatgaaac
tgcaatttattcatatcaggattatcaataccatatttttgaaaaagccgtttctgtaa
tgaaggagaaaaactcaccgaggcaggttccataggatggcaagatcctgggtatcggctg
cgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaataagg
ttatcaagtgagaaatcaccatgagtgcgactgaatccgggtgagaatggcaaaagt
atgcatttctttccagacttggtcaacaggccagccattacgctcgtcatcaaaatcac
tcgcatcaacaaaccgttattcattcgtgattgcgcctgagcgagacgaaatacgcga
tcgctgttaaaaggacaattacaaacaggaatcgaatgcaaccggcgcaggaacactgc
cagcgcatacaaatattttcacctgaatcaggatattcttctaataacctggaatgctg
ttttcccggggatcgcagtggtgagtaacctatgcatacaggagtaggataaaatgc
ttgatggtcggaagaggcataaaatccgctcagccagtttagtctgacctctcatctgt
aacatcattggcaacgctaccttgccatgtttcagaaacaactctggcgcatcgggt
tcccatacaatcgatagattgtcgacactgattgcccagacattatcgcgagcccattta
taccatataaatcagcatccatgttggaatttaatcgcggcctagagcaagacgtttc
ccgttgaatatggctcataacaccccttgattactgtttatgtaagcagacagtttta
ttgttcatgacaaaaatcccttaacgtgagttttcggtccactgagcgtcagacccgt
agaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgc
aaacaaaaaaaccaccgctaccagcgggtggtttgtttgcccggatcaagagctaccaact
ctttttccgaaggtaactggcttcagcagagcgcagataccaaatactgtccttctagt
gtagccgtagttaggccaccacttcaagaactctgtagcaccgcctacatacctcgctc
tgctaatactgttaccagtggtgctgccagtggcgataagtctgtgtttaccgggttg
gactcaagacgatagttaccggataaaggcgcagcggcggggtgaacgggggggttcgtg
cacacagcccagcttggagcgaacgcctacaccgaactgagatacctacagcgtgagc
tatgagaaagcgcacgccttccgaaggagaaaggcggacaggtatccggtaagcggc
aggggtcggaacaggagagcgcacaggggagcttccagggggaaacgcctgggtatcttta
tagtcctgtcgggtttcgccacctctgacttgagcgtcgatttttgatgctcgtcag
gggggaggagcctatggaaaaacgcagcaacgcggcctttttacggttcctggccttt
tgctggccttttgctcacatgttctttcctgcgttatcccctgattctgtggataaccg
tattaccgcctttgagtgcgtgataccgctcgccgcagccgaacgaccgagcgcagcg
agtcagtgcgcgaggaagcggaagagcgcctgatgcggtattttctccttacgcatctg
tgccgtattttcacaccgcataatgggtgcaactctcagtaaatctgctctgatgccgca
tagttaagccagtatatactccgctatcgctacgtgactgggtcatggctgcgccccga
caccgcgaacacccgctgacgcgcctgacgggcttgctctgctcccggcatccgctta
cagacaagctgtgaccgtctccgggagctgcatgtgtcagagggttttaccgctcatcac

Figure 81B

Nucleotide sequence of *P. tremuloides* pET24a

cgaaacgcgcgagggcagctgcggtaagctcatcagcgtggtcgtgaagcgattcacag
atgtctgectgttcatccgcgtccagctcgttgagtttctccagaagcgttaatgtctg
gcttctgataaagcgggcatgttaagggcgggtttttcctgtttggtcactgatgcct
ccgtgtaagggggatttctgttcatgggggtaatgataccgatgaaacgagagaggatg
ctcacgatacgggttactgatgatgaacatgcccggttactggaacgttgtgagggtaa
acaactggcgggtatggatgcggcgggaccagagaaaaatcactcagggtcaatgccagc
gcttcgttaatacagatgtagggtgtccacagggtagccagcagcatcctgcgatgcag
atccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacacg
gaaaccgaagaccattcatgttggtgctcaggtcgcagacgttttgacagcagcgtcgc
ttcacgttcgctcgcgtatcgggtgattcattctgctaaccagtaaggcaaccccgccag
cctagccgggtcctcaacgacaggagcagcatcatgcgcacccgtggggccgcatgcc
ggcgataatggcctgcttctcgcggaaacgtttggtggcgggaccagtgcgaaggctt
gagcgagggcggtgcaagattccgaataccgcaagcgacaggccgatcatcgtcgcgtc
cagcgaaagcgggtcctcgcggaaaatgaccagagcgtgccggcacctgtcctacgag
ttgcatgataaagaagacagtcataagtgcggcgacgatagtcatgccccgcgcccacc
ggaaggagctgactgggttgaaggctctcaaggcatcggtcgagatcccggtgcctaa
tgagtgagctaaactacattaattgcgttgcgctcactgcccgctttccagtcgggaaa
cctgtcgtgccagctgcattaatgaatcggccaacgcgcggggagaggcgggttgcgta
ttgggcgcaggggtggtttttcttttcaccagtgagacgggcaacagctgattgccctt
caccgcctggccctgagagagttgcagcaagcgggtccacgctggtttgccccagcaggc
gaaaatcctgtttgatggtggttaacggcgggatataacatgagctgtcttcggtatcg
tcgtatcccaactaccgagatataccgcaccaacgcgcagcccgactcggtaatggcgcg
cattgcgcccagcgccatctgatcgttggaaccagcatcgagtgggaaacgatgccct
cattcagcattttgatggtttggtgaaaaccggacatggcactccagtcgccttcccg
tccgctatcggtgaatttgattgcgagtgagatatttatgccagccagccagacgcag
acgcgcgagacagaacttaatgggcccgcctaacagcgcgatttgctggtgacccaatg
cgaccagatgctccacgcccagtcgcgtaccgtcttcatgggagaaaataatactgttg
atgggtgtctggtcagagacatcaagaaataacgcgggaacattagtcaggcagcttc
cacagcaatggcatcctggtcatccagcggatagttaatgatcagcccactgacgcgtt
gcgcgagaagattgtgcaccgcccgtttacaggcttcgacgcccgttcgttctaccatc
gacaccaccacgctggcaccagttgatcggcgcgagatttaatcgccgcgacaatttg
cgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacgactgtttgc
ccgccagttgttggtgccacgcggttggaatgtaattcagctccgccatcgccgcttcc
actttttcccgcttttcgcagaaacgtggctggcctggttcaccacgcgggaaacggt
ctgataagagacaccggcataactctgcgacatcgtataacgttactggtttcacattca
ccacctgaattgactctcttccgggcgctatcatgccataccgcgaaagggtttgcg
cattcgatggtgtccgggatctcgacgctctcccttatgcgactcctgcattaggaagc
agcccagtagtaggttgaggccgttgagcaccgcccgcgcaagggaatggtgatgcaag
gagatggcgcccaacagtcctccggccacggggcctgccaccatacccacgcgcgaaaca
agcgtcatgagcccgaagtggcgagcccgatcttcccatcggtgatgtcggcgatat
aggcgccagcaaccgcacctgtggcgccgggtgatgccggccacgatgcgtccggcgtag
aggatcgagatctcgatcccgcgaaattaatacgaactcactatagggggaattgtgagcg
gataacaattcccctctagaaataattttgtttaactttaagaaggagatatacatatg

Figure 81C

Nucleotide sequence of *P. tremuloides* pET24a

cgttgtagcgtgtccaccgaaaatgtgtctttctctgaaactgaaaccgaaacgcgtcg
ttctgcgaactacgaacctaacagctgggactatgattacctgctgtcctccgacacgg
acgagtcctatcgaagtacacaaagacaaagcgaaaaagctggaagccgaagttcgtcgc
gagattaataacgaaaaagcagaattttctgacctgctggaactgattgacaacgtcca
gcgcctgggcctgggttaccgtttcgagtctgatatccgtcgtgcgctggatcgcttcg
tttccctccggcggttcgatggcgtaaccaagacttccctgcacggtacggcactgtct
ttccgtctgctgcgtcaacacgggttttgaggtttctcaggaagcgttcagcggcttcaa
agacaaaacggcaacttccctggagaacctgaaggaagatatcaaagctatcctgagcc
tgtacgaggccagcttccctggctctggaaggcgaaaacatcctggacgaggcgaaaggtt
ttcgcaatctctcatctgaaagaactgtctgaagaaaagatcggtaaagagctggcaga
acaggtgtcccatgcaactggaactgccactgcatcgccgtactcagcgtctggaagcag
tatggtctatcgaggcctaccgtaaaaaggaggacgcgaaccaggttctgctggagctg
gcaattctggattacaacatgatccagtctgtataccagcgtgatctgcgtgaaacgtc
ccgttgggtggcgtcgtgtgggtctggcgaccaaactgcactttgctcgtgaccgcctga
ttgagagcttctactgggcgtgggtgtagcattcgaaccgcaataactccgactgccgt
aactccgtcgcaaaaatgttttcttctgtaaccattatcgacgatatctacgatgtata
cggcacccctggacgaactggagctgtttactgatgcagttgagcgttgggacgtaaacg
ccatcaacgacctgccggattacatgaaactgtgctttctggctctgtataaacactatt
aacgaaatcgctacgacaacctgaaagataaaggtgagaacatcctgccgtatctgac
caaagcctgggtgacctgtgcaacgctttcctgcaagaagccaagtggctgtacaaca
aatctactccgacctttgacgactacttcggcaacgcgatggaaatcctcttctggcccc
ctgcaactgatcttcgcttacttcgctgtcgtgcagaacattaaaaaggaagagatcga
aaacctgcaaaaataccatgacatcatctctcgtccttcccatatcttccgtctgtgca
atgacctggctagcgcgtctgcggaaattgcgcgtgggtgaaaccgcaaatagcgtttct
tggtacatgcgcaactaaaggtatctccgaagaactggctaccgaaagcgtgatgaatct
gatcgatgaaacctggaaaaagatgaacaaggaaaaactgggtggtagcctgttcgcga
aacggttcgtggaaaccgcgatcaacctggcacgtcaatctcactgcacttatcataac
ggcgacgcgcataacctctccgatgagctgaccgcgaaacgcgttctgtctgtaatcac
tgaaccgattctgccgtttgaacgctaaggatccgaattcgagctccgtcgacaagctt
ggggccgcactcgagcaccaccaccaccactgagatccggctgctaacaaagcccg
aaaggaagctgagttggctgctgccaccgctgagcaataactagcataaaccccttgggg
cctctaaacgggtcttgaggggttttttgcgtgaaaggaggaactatatccggat

Figure 82

Amino acid sequence of *P. trichocharpa* IspS (SEQ ID NO:160)

MRCSVSTENVSFTEETETETRRSANYEPNSWDYDYLLSSDTDESIEVYKDKAKKLEAEVR
REINNEKAEFLTLLELIDNVQRLGLGYRFESDIRRALDRFVSSGGFDAVTKTSLHATAL
SFRLLRQHGFVVSQEAQSGFKDQNGNFLENLKEDIKAILSLYEASFLALEGENILDEAK
VFAISHLKELSEEKIGKDLAEQVNHAEPLHRRRTQRLEAVLSIEAYRKEDADQVLE
LAILDYNNMIQSVYQRDLRETSRWRRRVGLATKLHFARDRLIESFYWAVGVAFEPQYSDC
RNSVAKMFSFVTIIDDYDVYGTLDLELFTNAVERWDVNAIDDLDPDYMKLCFLALYNT
INEIAYDNLKEKGENILPYLTAWADLCNAFLQEAQWLYNKSTPTFDEYFGNAWKSSSG
PLQLVFAYFAVVQNIKKEEIENLQKYHDIISRP SHIFRLCNDLASASAEIARGETANSV
SCYMRKGI SEELATESVMNLI DETWKKMNKEKLGGS LFAKPFVETAINLARQSHCTYH
NGDAHTSPDELTRKRVLSVITEPILPFER

Figure 83A

Nucleotide sequence of *P. trichocharpa* pET24a (SEQ ID NO:161)

tggcgaatgggacgcgcacctgtagcggcgcatctaagcgcggcggtgtggtggttacgc
gcagcgtgaccgctacacttgccagcgccctagcgcgcgctcctttcgctttcttccct
tcctttctcgccacgttcgcccggctttccccgtcaagctctaaatcggggggtcccttt
agggttccgatttagtgctttacggcacctcgacccccaaaaacttgattagggatgatg
gttcacgtagtgggcatcgccctgatagacgggttttcgccccttgacgttgagtc
acgttctttaatagtggactcttggtccaaactggaacaactcaaccctatctcgg
ctattcttttgatttataaagggttttgccgatttcggcctattgggttaaaaaatgagc
tgatttaacaaaaatttaacgcgaattttaacaaaatattaacgtttacaatttcagg
ggcacttttcggggaaatgtgcgcggaacccctatttggttatttttctaaatacattc
aaatatgtatccgctcatgaattaatcttagaaaaactcatcgagcatcaaatgaaac
tgcaatttattcatatcaggattatcaataccatattttgaaaaagccgtttctgtaa
tgaaggagaaaaactcacgagggcaggtccataggatggcaagatcctggatcggctc
cgattccgactcgtccaacatcaatacaacctattaatttcccctcgtcaaaaataagg
ttatcaagtgagaaatcacatgagtgacgactgaatccgggtgagaatggcaaaagttt
atgcatttctttccagacttggttaacaggccagccattacgctcgtcatcaaatcac
tcgcatcaaccaaaccgttattcattcgtgattgcgcctgagcgagacgaaatacgcga
tcgctgttaaaaggacaattacaaacaggaatcgaatgcaaccggcgaggaacactgc
cagcgcacacaaatattttcacctgaatcaggatattcttctaatacctggaatgctg
ttttcccggggatcgagtggtgagtaacctgcatcatcaggagtacggataaaatgc
ttgatggtcggaagaggcataaattccgtcagccagtttagtctgaccatctcatctgt
aacatcattggcaacgctacctttgccatgtttcagaaacaactctggcgcatcgggct
tcccatacaatcgatagattgtcgcacctgattgcccagacattatcgcgagcccattta
taccatataaatcagcatccatgttggaatttaatcgcggcctagagcaagacgtttc
cgttggaatattggtcataacaccccttgattactgtttatgtaagcagacagtttta
ttgttcatgacaaaaatcccttaacgtgagttttcgttccactgagcgtcagacccgt
agaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgc
aaacaaaaaaaccacgctaccagcgggtggtttgtttgccggatcaagagctaccaact
ctttttccgaaggtaactggcttcagcagagcgcagataccaaatactgtccttctagt
gtagcgttagttaggccaccacttcaagaactctgtagcaccgcctacatacctcgtc
tgctaactcctgttaccagtggctgctgccagtggcgataagtctgtcttaccgggttg
gactcaagacgatagttaccggataaggcgcagcggtcgggctgaacgggggggttcgtg
cacacagcccagcttgagcgaacgacctacaccgaactgagatacctacagcgtgagc
tatgagaaaagcgccacgcttccgaaggagaaaaggcgacaggtatccggttaagcggc
agggtcggaacaggagagcgcacgagggagcttccagggggaaacgcctggtatcttta
tagtcctgtcgggtttcgccacctctgacttgagcgtcgatttttgatgctcgtcag
ggggggcggagcctatggaaaaacgccagcaacgcggcctttttacggttccctggccttt
tgctggccttttgctcacatgttctttcctgcgttatcccctgattctgtggataaccg
tattaccgcctttgagtgagctgataccgctcgcgcagccgaacgaccgagcgcagcg
agtcagtgagcaggaagcggaagagcgctgatgcggtattttctccttacgcatctg
tgcggtatttcacaccgcatataggtgcaactctcagtacaatctgctctgatgccga
tagttaagccagtatacactccgctatcgctacgtgactgggtcatggctgcgccccga
caccgcgaacacccgctgacgcgcctgacgggcttgctgctcccggcatccgctta
cagacaagctgtgaccgtctccggagctgcatgtgtcagagggttttaccgctcatcac

Figure 83B

Nucleotide sequence of *P. trichocharpa* pET24a

cgaaacgcgcgagggcagctgcggtaaagctcatcagcgtggtcgtgaagcgattcacag
atgtctgcctgttcatccgcgtccagctcggttgagtttctccagaagcggtaatgtctg
gcttctgataaagcgggccatgttaagggcggttttttctgtttggtcactgatgcct
ccgtgtaagggggatttctgttcatgggggtaatgataccgatgaaacgagagaggatg
ctcacgatacgggttactgatgatgaacatgcccggttactggaacgttgtgagggtaa
acaactggcggtatggatgcggcgggaccagagaaaaatcactcagggccaatgccagc
gcttcgttaatacagatgtagggtgttccacagggtagccagcagcatcctgcgatgcag
atccggaacataatggtgcagggcgctgacttccgcgtttccagactttacgaaacacg
gaaaccgaagaccattcâtgttggtgctcaggtcgagacggttttgagcagcagtcgc
ttcacgttcgctcgcgtatcgggtgattcattctgctaaccagtaaggcaaccccgcag
cctagccgggtcctcaacgacaggagcagcatcatgcgcacccgtggggccgcatgcc
ggcgataatggcctgcttctcgcgaaacgtttggtggcgggaccagtgcgaaggctt
gagcgagggcggtgcaagattccgaataccgcaagcgacaggccgatcatcgtcgcgtc
cagcgaaagcggtcctcgcgaaatgacccagagcgctgccggcacctgtcctacgag
ttgcatgataaagaagacagtcataagtgcggcgacgatagtcatgccccgcgcccacc
ggaaggagctgactgggttgaaggctctcaagggcatcggtcgagatccgggtgcctaa
tgagtgagctaaactacattaattgcgttgcgctcactgcccgctttccagtcgggaaa
cctgtcgtgccagctgcattaatgaatcggccaaacgcgcggggagaggcggtttgcgta
ttgggcgccagggtggtttttcttttaccagtgagacgggcaacagctgattgccctt
caccgcctggccctgagagagtgcagcaagcggtccacgctggtttgccccagcaggc
gaaaaatcctgtttgatgggtggttaacggcgggatataacatgagctgtcttcggatcg
tcgtatcccactaccgagatatccgcaccaacgcgcagcccgactcggtaatggcgcg
cattgcgcccagcgccatctgatcgttggaaccagcatcgagtgggaaacgatgccct
cattcagcatttgcatggtttgttgaaaaccggacatggcactccagtcgccttcccgt
tccgctatcggtgaatttgattgcgagtgagatatttatgccagccagccagacgcag
acgcgccgagacagaactaatgggcccgcataacagecgcatttgctggtgacccaatg
cgaccagatgctccacgcccagtcgcgtaccgtcttcatgggagaaaaataactgttg
atgggtgtctggtcagagacatcaagaaataacgcgggaacattagtgcaggcagcttc
cacagcaatggcatcctgggtcatccagcggatagttaatgatcagcccactgacgcgtt
gcgcgagaagattgtgcaccgcgcgtttacaggcttcgacgcgcgttcgttctaccatc
gacaccaccacgctggcaccagttgatcggcgcgagatttaatcgccgcgacaatttg
cgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacgactgtttgc
ccgccagttgttggtgccacgcggttggaatgtaattcagctccgccatcgccgcttcc
actttttcccgcgttttcgcagaaacgtgggtggcctggttcaccacgcgggaaacggt
ctgataagagacaccggcataactctgcgacatcgtataacgttactggtttcacattca
ccaccctgaattgactctcttccgggcgctatcatgccataaccgcgaaagggttttgcg
cattcgatggtgtccgggatctcgacgctctcccttatgcgactcctgcattaggaagc
agcccagtagtaggttgaggccgttgagcaccgcgcgcgcaaggaatggtgcatgcaag
gagatggcgcccaacagtcctccggccacggggcctgccaccatacccacgcgcgaaaca
agcgcctcatgagcccgaagtggcgagcccgatcttccccatcggatgctcggcgatat
aggcgccagcaaccgcacctgtggcgccggtgatgccggccacgatgcgtccggcgtag
aggatcgagatctcgatccgcgaaatttaatacagactcactatagggaattgtgagcg
gataacaattcccctctagaaataattttgtttaactttaagaaggagatatacatatg

Figure 83C

Nucleotide sequence of *P. trichocharpa* pET24a

catatgcgttgtagcgtgtccaccgaaaatgtgtctttcaccgaaactgaaaccgaaac
gcgtcggttctgcgaactacgaacctaacagctgggactatgattacctgctgcctccg
acacggacgagtcctatcgaagtatacaaagacaaagcgaaaaagctggaagccgaagt
cgtcgcgagattataaacgaaaaagcagaatttctgacctgctggaactgattgacaa
cgccagcgccctgggctgggttacggttctcgagtctgatccgctcgctgcgtggatc
gcttcgtttcctccggcggcctcgatgcggttaaccaagacttccctgcacgcgacggca
ctgtctttccgtctgctgcgtcaacacgggttttgaggtttctcaggaagcgttcagcgg
cttcaaagacaaaaacggcaacttccctggagaacctgaaggaagatatcaaagctatcc
tgagcctgtacgagggccagcttccctggctctggaaggcgaaaacatccctggacgaggcg
aaggttttcgcaatctctcatctgaaagaactgtctgaagaaaagatcggtaaagatct
ggcagaacaggtgaacctgcaactggaactgccactgcatcgccgtactcagcgtctgg
aagcagtactgtctatcgaggcctaccgtaaaaaggaggacgcggatcaggttctgctg
gagctggcaattctggattacaacatgatccagtctgtataccagcgtgatctgcgtga
aacgtcccggttgggtggcgtcgtgtgggtctggcgaccaaactgcactttgctcgtgacc
gcctgattgagagcttctactgggcgctgggtgtagcattcgaaccgcaatactccgac
tgccgtaaactccgtcgcaaaaatgttttcttctgtaaccattatcgacgatctacga
tgtatacggcacccctggacgaactggagctgttactaacgcagttgagcgttgggacg
taaacgccatcgacgatctgcgggattacatgaaactgtgctttctggctctgataac
actattaacgaaatcgctacgacaacctgaaagaaaaaggtgagaacatcctgccgta
tctgaccaaagcctgggtgacctgtgcaacgctttcctgcaagaagccaagtggctgt
acaacaaatctactccgacctttgacgaatacttcggcaacgcattgaaatcctcttct
ggcccgctgcaactggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaaga
gatcgaaaacctgcaaaaataccatgacatcatctctcgtccttcccatatcttcgctc
tgtgcaatgacctggctagcgcgtctgcggaaattgcgcgtggtgaaaccgcaaatagc
gtttcttgttacatgcgcactaaaggtatctccgaagaactggctaccgaaagcgtgat
gaatctgatcgatgaaacctggaaaaagatgaacaaggaaaaactgggtggtagcctgt
tcgcgaaaccggttcgtggaaaccgcgatcaacctggcacgtcaatctcactgcacttat
cataacggcgacgcgcataacctctccggatgagctgacccgcaaaccggttctgtctgt
aatcactgaaccgattctgccgtttgaacgctaaggatccgaattcgagctccgtcgac
aagcttgccggccgcactcgagcaccaccaccaccactgagatccggctgctaacaa
agcccgaaaggaagctgagttggctgctgccaccgctgagcaataactagcataacccc
ttggggcctctaaacgggtcttgaggggttttttgctgaaaggaggaactatatccgga
t

Figure 84A




Figure 84B




Figure 84C




Figure 85

Amino acid sequence of IspS variant P. albaTRC-pET200 (SEQ ID NO:162)

MRGSHHHHHGMA SMTGGQQMGRDLYDDDDKDH PFTMRRSANYEPNSWDYDYLLSSD TD
ESIEVYKDKAKKLEAEVRREINNEKAEFLT LLELIDNVQRLGLGYRFESDIRGALDRFV
SSGGFDAVTKTSLHGTALSFRLLRQHGF EVSQEAFSGFKDQNGNFLENLKEDIKAILSL
YEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVNHALELPLHRR TQRLEAV
WSIEAYRKKEDANQV LLELA ILDYNMIQSVYQRDLRETSRWRRVGLATKLHFARDRLI
ESFYWAVGVAFEPQYSDCRNSVAKMFSFVTI IDDIYDVYGTLDLELFTDAVERWDVNA
INDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNK
STPTFDDYFGNAWKSSSGPLQLVFAYFAVVQNIKKEE IENLQKYHDTISRPSHIFRLCN
DLASASAEIARGETANSVSCYMR TKGISEELATESVMNLIDETWKKMNKEKLGGSLFAK
PFVETAINLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 86A

Nucleotide Sequence of pDu30 (SEQ ID NO:163)

aagggcgagctcaacgatccggctgctaacaaagcccgaaaggaagctgagttggctgc
tgccaccgctgagcaataactagcataaacccttggggcctctaaacgggtcttgagga
gttttttgcgtgaaaggaggaactatatccggatatcccgcaagaggcccggcagtagcg
gcataaccaagcctatgcctacagcatccagggtgacgggtgccgaggatgacgatgagc
gcattgttagatttcatcacgggtgcctgactgcgttagcaatttaactgtgataaact
accgcattaaagccttatcgatgataagctgtcaaactgagaatttaattcttgaagacg
aaagggcctcgtgatacgcctatttttataggttaatgtcatgataataatgggttctt
agacgtcagggtggcacttttcggggaaatgtgcgcggaaccctatttgtttattttc
taaatacattcaaataatgtatccgctcatgagacaataaccctgataaatgcttcaata
atattgaaaaaggaagagtatgattgaacaagatggattgcacgcaggttctccggccg
cttgggtggagaggctattcggctatgactgggcacaactgacaatcggctgctctgat
gccgcgctgttccggctgtcagcgcagggggcgcccggttcttttgcagaccgacct
gtccggtgccctgaatgaactgcaggacgaggcagcgcggctatcgtggctggccacga
cgggcgttccctgcgcagctgtgctcgacgttgtcactgaagcgggaagggactggctg
ctattgggcgaagtgcggggcaggatctcctgtcatctcaccttgcctcctgccgagaa
agtatccatcatggctgatgcaatgcggcggtgcatacgcctgatccggctacctgcc
cattcgaccaccaagcgaaacatcgcatcgagcgggcacgtactcggatggaagccggt
cttgtcgatcaggatgatctggacgaagagcatcaggggtcgcgccagccgaactgtt
cgccaggctcaaggcgcgcgatgcccgacggcgaggatctcgtcgtgacacatggcgatg
cctgcttggcgaatatcatggtggaaaatggccgcttttctggattcatcgactgtggc
cggtgggtgtggcggaaccgctatcaggacatagcgttggctaccgctgatattgctga
agagcttggcggcgaatgggtgaccgcttccctcgtgctttacgggtatcgccgctcccg
attcgcagcgcacgccttctatcgcccttcttgacgagttcttctgagcgggactctgg
gggtcgaaatgaccgaccaagcgacgcctaaactgtcagaccaagtttactcatatatac
tttagattgattttaaacttcattttttaattttaaaggatctaggtgaagatccttttt
gataatctcatgacaaaaatcccttaacgtgagttttcgttccactgagcgtcagaccc
cgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgct
tgcaaacaaaaaaaccacgctaccagcgggtggtttgtttgccggatcaagagctacca
actctttttccgaaggtaactggcttcagcagagcgcagataccaaatactgtccttct
agtgtagccgtagtttaggccaccacttcaagaactctgtagcaccgcctacatacctcg
ctctgctaatacctgttaccagtggctgctgccagtggcgataagtctgtcttaccggg
ttggactcaagacgatagttaccggataaaggcgcagcggctcgggctgaacggggggttc
gtgcacacagcccagcttggagcgaacgacctacaccgaactgagatacctacagcgtg
agctatgagaaagcgccacgcttcccgaaggagaaaggcggacaggtatccggttaagc
ggcaggggtcggaacaggagagcgcacgagggagcttccagggggaaacgcctggtatct
ttatagtccgtcgggtttccgccacctctgacttgagcgtcgatttttgtgatgctcgt
cagggggggcgagcctatggaaaaacgccagcaacgcggcctttttacgggttccctggcc
ttttgctggccttttgcctcacatgttcttctcctgcgttatcccctgattctgtggataa
ccgtattaccgcctttgagtgagctgataccgctcgccgcagccgaacgaccgagcgca
gcgagtcagtgagcgaggaagcggaagagcgcctgatgcggtattttctccttacgcat
ctgtgcggtatttcacaccgcaatggtgcaactctcagtacaatctgctctgatgccgca
tagttaagccagtatatactccgctatcgctacgtgactgggtcatggctgcgccccga
caccgcgaacacccgctgacgcgcctgacgggcttgtctgctcccggcatccgctta

Figure 86B

Nucleotide Sequence of pDu30

cagacaagctgtgaccgtctccgggagctgcatgtgtcagagggttttaccggtcatcac
cgaaacgcgcgagggcagctgcggtaaaagctcatcagcgtggtcgtgaagcgattcacag
atgtctgcctgttcatccgcgtccagctcgttgagtttctccagaagcgttaatgtctg
gcttctgataaagcgggccaatgttaagggcgggttttttctggttggtcactgatgcct
ccgtgtaagggggatttctgttcatgggggtaatgataccgatgaaacgagagaggatg
ctcacgatacgggttactgatgatgaacatgccgggttactggaacgttgtgagggtaa
acaactggcggatggatgcggcgggaccagagaaaaatcactcagggtcaatgccagc
gcttcgttaatacagatgtagggtgttccacagggtagccagcagcatcctgcgatgcag
atccggaacataatggtgcagggcgctgacttccgcgttccagactttacgaaacacg
gaaaccgaagaccattcatgttggtgctcaggtcgcagacgttttgcagcagcagtcgc
ttcacgttcgctcgcgtatcgggtgatcattctgctaaccagtaaggcaaccccgccag
cctagccgggtcctcaacgacaggagcagcatcatgcgcaccggtggccaggaccaac
gctgcccagatgcgcgcgtgcgggtgctggagatggcggacgcgatggatgttct
gccaagggttggtttgcgcattcacagttctccgcaagaattgattggctccaattctt
ggagtggatgaatccgttagcaggtgcccgcggcttccattcagggtcgagggtggcccg
ctccatgcaccgcgacgcaacgcggggagggcagacaaggatagggcgggcgccataat
ccatgccaacccgttccatgtgctgcgcgagggcgccataaatcccggtgacgatcagcg
gtccaatgatcgaagttagggtggttaagagccgcgagcgatccttgaagctgtccctga
tggtcgtcatctacctgcctggacagcatggcctgcaacgcggggcatccgatgcgcgc
ggaagcgagaagaatcataatggggaaggccatccagcctcgctcgcgaacgccagca
agacgtagcccagcgcgtcgcccgccatgccggcgataatggcctgcttctcgccgaaa
cggttggtggcgggaccagtgacgaaggcttgagcgaggcggtgcaagattccgaatac
cgcaagcgacaggccgatcatcgtcgcgtccagcgaaagcggctcctcgccgaaaatga
cccagagcgctgccggcacctgtcctacgagttgcatgataaagaagacagtcataagt
gcggcgacgatagtcatgccccgcgcccaccggaaggagctgactgggttgaaggctct
caagggcatcggtcgagatcccggtgcctaatgagtgaagtaacttacattaattgcgt
tgcgctcactgcccgtttccagtcgggaaacctgtcgtgccagctgcattaatgaatc
ggccaacgcgcggggagagggcggtttgcggtattgggcgccagggtgggttttcttttca
ccagtgcagcgggcaacagctgattgcccttcaccgcctggccctgagagagttgcagc
aagcgggtccacgctggtttgccccagcaggcgaaaatcctgttgatgggtggttaacgg
cgggatataacatgagctgtcttcggtatcgtcgtatcccactaccgagatatccgcac
caacgcgcagcccggactcggtaatggcgcgcattgcgcccagcgccatctgatcgttg
gcaaccagcatcgcagtggaacgatgccctcattcagcatttgcatgggttggtgaaa
accggacatggcactccagtcgccttcccggttcgcgtatcggtgaatttgattgcgag
tgagatatttatgccagccagccagacgcgagcgcgagacagaacttaatgggccc
gctaacagcgcgatttgctgggtgacccaatgcgaccagatgctccacgcccagtcgcgt
accgtcttcatgggagaaaaataatactgttgatgggtgtctggtcagagacatcaagaa
ataacgcgggaacattagtgcaggcagcttccacagcaatggcatcctggtcatccagc
ggatagttaatgatcagcccactgacgcgttgcgcgagaagattgtgcaccgcgcgttt
acaggttcgacgcgcgttctgttctaccatcgacaccaccacgctggcaccagttgat
cggcgcgagatttaatcgccgcgacaatttgcgacggcgcggtgcagggccagactggag
gtggcaacgccaatcagcaacgactgtttgcccgcagttgttggtgccacgcggttggg
aatgtaattcagctccgccatcgccgcttccacttttcccgcggttttcgcagaaacgt

Figure 86C

Nucleotide Sequence of pDu30

ggctggcctggttcaccacgcgggaaacggtctgataagagacacccgcataactctgcg
acatcgtataacggttactggtttcacattcaccaccctgaattgactctctccgggcg
ctatcatgccataccgcgaaaggttttgcgccattcgatgggtgtccgggatctcgacgc
tctcccttatgcgactcctgcattaggaagcagcccagtagtaggttgaggccgttgag
caccgcgcgcgaaggaatggtgcatgcaaggagatggcgcccaacagtccccggcca
cggggcctgccaccatacccacgcgcgaaacaagcgctcatgagcccgaagtggcgagcc
cgatcttccccatcgggtgatgtcggcgatataggcgccagcaaccgcacctgtggcgcc
ggtgatgccggccacgatgcgtcggcgtagaggatcgagatctcgatcccgcgaaatt
aatacgactcactataggggaattgtgagcggataaacaattcccctctagaaataattt
tgtttaactttaagaaggagatatacatatgcggggttctcatcatcatcatcatg
gtatggctagcatgactggtggacagcaaattgggtcgggatctgtacgacgatgacgat
aaggatcatcccttcaccatgcgtcgttctgcgaactacgaacctaacagctgggacta
tgattacctgctgtcctccgacacggacgagtcctcgaagtatacaaagacaaagcga
aaaagctggaagccgaagttcgtcgcgagattaataacgaaaaagcagaatttctgacc
ctgctggaactgattgacaacgtccagcgcctgggcctgggttacggtttcgagctctga
tatccgtggtgcgctggatcgcttcgtttcctccggcggcttcgatgcggtaaccaaga
cttccctgcacggtagcgcaactgtcttccgtctgctgcgtcaacacggttttgaggtt
tctcaggaagcgttcagcggcttcaaagacaaaacggcaacttccctggagaacctgaa
ggaagatatcaaagctatcctgagcctgtacgaggccagcttccctggctctggaaggcg
aaaacatcctggacgagggcgaaggttttcgcaatctctcatctgaaagaactgtctgaa
gaaaagatcggtaaaagagctggcagaacaggtgaacctgcaactggaactgccactgca
tcgccgtactcagcgtctggaagcagtatggtctatcgaggcctaccgtaaaaaggagg
acgcgaatcaggttctgctggagctggcaattctggattacaacatgatccagtctgta
taccagcgtgatctgcgtgaaacgtcccggttggtggcgtcgtgtgggtctggcgaccaa
actgcactttgctcgtgacgcctgattgagagcttctactgggcgtgggtgtagcat
tcgaaccgcaatactccgactgcgtaactcgtcgcaaaaatgttttcttctgtaacc
attatcgacgatatactacgatgtatacggcacccctggacgaactggagctgtttactga
tgcagttgagcgttgggacgtaaacgccatcaacgacctgccggattacatgaaactgt
gctttctggctctgtataacactattaacgaaatcgctacgacaacctgaaagataaa
gggtgagaacatcctgccgtatctgaccaaagcctgggctgacctgtgcaacgctttcct
gcaagaagccaagtggctgtacaacaaatctactccgaccttgacgactacttcggca
acgcattggaaatcctcttctggcccgtgcaactggtgttcgcttacttcgctgtcgtg
cagaacattaaaaaggaagagatcgaaaacctgcaaaaataccatgacaccatctctcg
tccttcccatatcttccgtctgtgcaatgacctggctagcgcgtctgcggaaattgcgc
gtggtgaaaccgcaaatagcgtttcttggtacatgcgcactaaaggatctccgaagaa
ctggctaccgaaagcgtgatgaatctgatcgatgaaacctggaaaaagatgaacaagga
aaaactgggtggtagcctgttcgcgaaaccgttcgtggaaccgcgatcaacctggcac
gtcaatctcactgcacttatcataacggcgacgcgcataacctctccggatgagctgacc
cgcaaacgcgttctgtctgtaatactgaaccgattctgccgtttgaacgctaa

Figure 87

Amino acid sequence of IspS variant *P.trem*TRC-pET200 (SEQ ID NO:164)

MRGSHHHHHGMASTGGQQMGRDLYDDDDKDHFFTMRRSANYEPNSWDYDYLLSSD
ESIEVHKDKAKKLEAEVRREINNEKAEFLTLELIDNVQRLGLGYRFESDIRRALDRFV
SSGGFDGVTKTSLHGTALSFRLLRQHGFVSQEAQSGFKDQNGNFLENLKEDIKAILSL
YEASFLALEGENILDEAKVFAISHLKELSEEKIGKELAEQVSHALELPLHRRRTQRLEAV
WSIEAYRKKEDANQVLELAILDYNMIQSVYQRDLRETSRWRRVGLATKLHFARDRLI
ESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDYDVYGTLDLELFTDAVERWDVNA
INDLPDYMKLCFLALYNTINEIAYDNLKDKGENILPYLTKAWADLCNAFLQEAKWLYNK
STPTFDDYFGNAWKSSSGPLQLIFAYFAVVQNIKKEIENLQKYHDIISRPSHIFRLCN
DLASASAEIARGETANSVSCYMRTKGI SEELATESVMNLIDETWKKMNKEKLGGSLFAK
PFVETAINLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 88A

Nucleotide sequence of pDu31 (SEQ ID NO:165)

cgtcgttctgcgaactacgaacctaacagctgggactatgattacctgctgtcctccga
cacggacgagtcctcgaagtacacaaagacaaagcgaaaaagctggaagccgaagttc
gtcgcgagatttaataacgaaaaagcagaatttctgacctgctggaactgattgacaac
gtccagegcctgggctgggttaccgtttcgagctctgatccgtcgtgcgctggatcg
cttcgtttcctccggcggttcgatggcgtaaccaagacttcctgcacggtagggcac
tgtctttccgtctgctgcgtcaacacgggttttgaggtttctcaggaagcgttcagcggc
ttcaaagacccaaacggcaacttcctggagaacctgaaggaagatatcaaagctatcct
gagcctgtacgagggcagcttcctggctctggaaggcgaaaacatcctggacgaggcga
aggttttcgcaatctctcatctgaaagaactgtctgaagaaaagatcggtaaaagagctg
gcagaacaggtgtcccatgcaactggaactgccactgcatcgccgtactcagcgtctgga
agcagtatggtctatcgaggcctaccgtaaaaaggaggacgcgaaccaggttctgctgg
agctggcaattctggattacaacatgatccagtctgtataccagcgtgatctgctgaa
acgtcccgttggtagcgtcgtgtgggtctggcgacccaaactgcactttgctcgtgaccg
cctgattgagagcttctactgggcgtgggtgtagcattcgaaccgcaatactccgact
gccgtaaactccgtcgcaaaaatgttttcttctgtaaccattatcgacgatatctacgat
gtatacggcaccctggacgaactggagctgttactgatgcagttgagcgttgggacgt
aaacgccatcaacgacctgccggattacatgaaactgtgctttctggctctgtataaca
ctattaacgaaatcgccctacgacaacctgaaagataaagggtgagaacatcctgccgtat
ctgacccaaagcctgggctgacctgtgcaacgctttcctgcaagaagccaagtggctgta
caacaaatctactccgacctttgacgactacttcggcaacgcattgaaatcctctctg
gcccgtgcaactgatcttcgcttacttcgctgtcgtgcagaacattaaaaaggaagag
atcgaaaacctgcaaaaataccatgacatcatctctcgtccttcccatatcttcgctct
gtgcaatgacctggctagcgcgtctgcggaaattgcgcgtggtgaaaccgcaaatagcg
tttcttggttacatgcgcactaaaggatctccgaagaactggctaccgaaagcgtgatg
aatctgatcgatgaaacctggaaaaagatgaacaaggaaaaactgggtggtagcctggt
cgcgaaaccgttcgtggaaaccgcgatcaacctggcacgtcaatctcactgcacttatac
ataacggcgacgcgcatacctctccggatgagctgacccgcaaacgcgttctgtctgta
atcactgaaccgattctgccgtttgaacgctaaaagggcgagctcaacgatccggctgc
taacaaagcccgaaggaagctgagttggctgctgccaccgctgagcaataactagcat
aacccttggggcctctaaacgggtcttgaggagtttttgctgaaaggaggaactata
tccggatatcccgaagaggcccgagcaccggcataaccaagcctatgcctacagca
tccaggggtgacgggtgccgaggatgacgatgagcgcattggttagatttcatacaggtgc
ctgactgcgttagcaatttaactgtgataaactaccgcattaaagcttatcgatgataa
gctgtcaaacatgagaatttaattcttgaagacgaaagggcctcgtgatagcctatctt
tataggttaattgtcatgataataatgggtttcttagacgtcaggtggcacttttcgggga
aatgtgcgcgggaacccctatctgtttatcttctaaatacattcaaataatgtatccgct
catgagacaataacccctgataaatgcttcaataatattgaaaaaggaagagtatgattg
aacaagatggattgcacgcaggttctccggccgcttgggtggagaggctattcggctat
gactgggcacaaactgacaatcggtctctgatgccgccgtgttcgggtgtcagcgca
ggggcgcccggttcttttctgcaagaccgacctgtccgggtgcctgaatgaactgcagg
acgaggcagcgcggctatcgtggctggccacgacgggcgttccttgcgcagctgtgctc
gacgttgtcactgaagcgggaagggaactggctgctattgggcgaagtgccggggcagga
tctcctgtcatctcaccttgctcctgccgagaaagtatccatcatggctgatgcaatgc

Figure 88B

Nucleotide sequence of pDu31

ggcggctgcatacgccttgatccggctacctgccattcgaccaccaagcgaaacatcgccagcgggcagcgtactcggatggaagccggctcttgatcgatcaggatgatctggacgaagagcatcaggggctcgcgccagccgaactgttcgccaggctcaaggcgcgcatgcccgacggcgaggatctcgtcgtgacacatggcgatgcctgcttgccgaatatcatggaggaaatggcgcgttttctggattcatcgactgtggccggctgggtgtggcggaccgctatcaggacatagcgttgggtaccgctgatattgctgaagagcttggcggcggaatgggctgaccgcttcctcgtgctttacgggtatcgccgctcccgatctgcagcgcatcgcttctatcgcttcttgacgagttcttctgagcgggactctgggggtcgaaatgaccgaccaagcgacgcctaactgtcagaccaagtttactcatatatactttagattgatttaaaacttcatttttaatttaaaaaggatctagggtgaagatcctttttgataatctcatgaccaaaatcctttaacgtgagttttcgttccactgagcgtcagaccccgtagaaaagatcaaaggatcttcttgagatccttttttctgcgcgtaatctgctgcttgcaacaaaaaaaccaccgctaccagcgggtggtttggttgccggatcaagagctaccaactccttttccgaaggtaactggcttcagcagagcgcagataccaaatactgtccttctagtgtagccgtagttaggccaccacttcaagaactctgtagcaccgcctacataacctcgctctgctaactcctgttaccagtggctgctgccagtggcgataagtctgtcttaccgggttggaactcaagacgatagttaccgga taaggcgagcgggtcgggctgaacgggggggttcgtgcacacagcccagcttggaagcgaa cgacctacaccgaactgagataacctacagcgtgagctatgagaaagcgccacgcttccc gaagggagaaaggcggacaggatccggtaagcggcagggtcggaacaggagagcgcacgagggagcttccagggggaaacgcctggtatctttatagtctgtcgggtttcgccacc tctgacttgagcgtcgatttttgatgctcgtcaggggggcggagcctatggaaaaac gccagcaacgcggccttttacgggttccctggccttttgctggccttttgctcacatgtt ctttctcgttatccctgattctgtggataaccgtattaccgcctttgagtgagctg ataccgctcgccgcagccgaacgaccgagcgcagcagtcagtgagcgaggaagcggaa gagcgctgatcgggtattttctccttacgcctctgtcgggtatttcacaccgcaatggtgcactctcagtacaatctgctctgatgccgcatagttaagccagtatacactccgcta tgcctacgtgactgggtcatggctgcgccccgacacccgcgaacacccgctgacgcgcc ctgacgggcttgctgctcccggcatccgcttacagacaagctgtgaccgtctccggga gctgcatgtgtcagagggttttaccgctcatcaccgaaacgcgcgagggcagctgcggtaa agctcatcagcgtggctcgtgaagcgattcacagatgtctgcctgttcatccgcgtccag ctggttagtttctccagaagcgttaatgtctggcttctgataaagcggggccatgttaa gggcgggtttttcctggttggtcactgatgcctccgtgtaaggggggatttctgttcatg ggggtaatgataccgatgaaacgagagaggatgctcacgatacgggttactgatgatga acatgcccggttactggaacgttgtgagggtaaacaaactggcgggtatggatgcggcggg accagagaaaaatcactcagggtcaatgccagcgttcgttaatacagatgtagggtgtt ccacagggtagccagcagcatcctgcgatgcagatccggaacataatggtgcaggggcgc tgacttccgcgtttccagactttacgaaacacggaaaccgaagaccattcatgttgttg ctcaggtcgagacgttttgagcagcagtcgcttcacgttcgctcgcgatcggatgat tcatctgctaaccagtaaggcaaccccgccagcctagccgggtcctcaacgacaggag cacgatcatgcgcaccgctggccaggaccaacgctgcccgagatgcgcgcgctgcggc tgcgtggagatggcggacgcgatggatattgttctgccaagggttggtttgcgcattcaca gtttctccgcaagaattgattgggtccaattcttgaggatgggtgaatccgttagcgagggtg ccgccggttccattcaggtcgaggtggccgggtccatgcaccgcgacgcaacgcggg

Figure 88C

Nucleotide sequence of pDu31

gaggcagacaaggtatagggcggcgcctacaatccatgccaacccgttccatgtgctcg
ccgagggcggcataaatcgccgtgacgatcagcgggtccaatgatcgaagttaggctggta
agagccgcgagcgatccttgaagctgtccctgatggtcgtcatctacctgcctggacag
catggcctgcaacgcgggcatcccgatgccgcgggaagcgagaagaatcataatgggga
aggccatccagcctcgctcgcgaaacgccaagacgtagcccagcgcgtcgggccgc
atgccggcgataatggcctgcttctcgccgaaacgtttgggtggcgggaccagtgacgaa
ggcttgagcgagggcggtgcaagattccgaataccgcaagcgacaggccgatcatcgtcg
cgctccagcgaaagcggtcctcgccgaaaatgaccagagcgctgccggcacctgtcct
acgagttgcatgataaagaagacagtcataagtgcggcgacgatagtcagccccgcgc
ccaccggaaggagctgactgggttgaaggctctcaagggcatcggtcgagatcccgggtg
cctaatgagtgagctaaacttacattaattgcgttgcgtcactgcccgtttccagtcg
ggaaacctgtcgtgccagctgcattaatgaatcgccaacgcgcggggagaggcggttt
gcgtattgggcgccagggtgggtttttcttttaccagtgagacgggcaacagctgattg
cccttcaccgcctggccctgagagagttgcagcaagcgggtccacgctggtttgccccag
caggcgaaaatcctgtttgatgggtgggttaacggcgggatataacatgagctgtcttcgg
tatcgtcgtatcccactaccgagatataccgcaccaacgcgcagcccggactcggtaatg
gcgcgcattgcgcccagcgccatctgatcgttggcaaccagcatcgagtggaacgat
gccctcattcagcatttgcatgggtttggtgaaaaccggacatggcactccagtcgcctt
cccgttccgctatcggtgaatttgattgcgagtgagatatttatgccagccagccaga
cgcagacgcgcgcgagacagaacttaatgggcccgttaacagcgcgatttgctggtgacc
caatgcgaccagatgctccacgcccagtcgcgtaccgtcttcatgggagaaaataatac
tggttgatgggtgtctggtcagagacatcaagaaataacgccggaacattagtgcaggca
gcttccacagcaatggcatcctggtcacccagcgatagttaatgatcagcccactgac
gcgttgcgcgagaagattgtgcaccgcgcgtttacaggcttcgacgcgcgttcgttcta
ccatcgacaccaccaacgctggcaccagttgatcggcgcgagatttaatcgccgcgaca
at ttgcgacggcgcggtgcagggccagactggaggtggcaacgcgcaatcagcaacgactg
tttgcccgcagttgttgtgccacgcggttgggaatgtaattcagctccgccatcgccg
cttccactttttcccgcttttcgcgaaacgtggctggcctggttcaccacgcgggaa
acggtctgataagagacacccgcatactctgcgacatcgtataacgttactggtttcac
attcaccacccctgaattgactctcttcggggcgctatcatgccataccgcgaaagggttt
tgcgccattcgatgggtgtccgggatctcgacgctctcccttatgcgactcctgcattag
gaagcagcccagtagtaggttgaggccgttgagcaccgcgcgcgcaagggaatggtgcat
gcaaggagatggcgcccaacagtcgcccgccacggggcctgccaccataccacgcgcg
aaacaagcgctcatgagccogaagtggcgagcccgatcttccccatcggtgatgtcggc
gatatagggcgccagcaaccgcacctgtggcgccggtgatgccggccacgatgcgtccgg
cgtagaggatcgagatctcgatcccgcgaaatttaatacgactcactataggggaattgt
gagcggataacaattcccctctagaaataattttgtttaactttaagaaggagatatac
atatcggggggttctcatcatcatcatcatcatggtagcatgactggtggacag
caaatgggtcgggatctgtacgacgatgacgataaggatcatcccttcaccatg

Figure 89

Amino acid sequence of IspS variant P.trichTRC-pET200 (SEQ ID NO:166)

MRGSHHHHHHGMASMTGGQQMGRDLYDDDDKDHFFTMRRSANYEPNSWDYDYLLSSDTD
ESIEVYKDKAKKLEAEVRREINNEKAEFLTLELIDNVQRLGLGYRFESDIRRALDRFV
SSGGFDAVTKTSLHATLSFRLLRQHGFVSQEAFSGFKDQNGNFLENLKEDIKAILSL
YEASFLALEGENILDEAKVFAISHLKELSEEKIGKDLAEQVNHAEPLHRRRTQRLEAV
LSIEAYRKKEADQVLLELAILDYNMIQSVYQRDLRETSRWRRVGLATKLHFARDRLI
ESFYWAVGVAFEPQYSDCRNSVAKMFSFVTIIDDIDYDVYGTLDLELFTNAVERWDVNA
IDDLPDYMKLCFLALYNTINEIAYDNLKEKGENILPYLTKAWADLCNAFLQEAKWLYNK
STPTFDEYFGNAWKSSSGPLQLVFAYFAVVQNIKKEIENLQKYHDIISRPSHIFRLCN
DLASASAEIARGETANSVSCYMRTKGISEELATESVMNLIDETWKKMNKEKLGGSLFAK
PFVETAINLARQSHCTYHNGDAHTSPDELTRKRVLSVITEPILPFER

Figure 90A

Nucleotide Sequence of pDu32 (SEQ ID NO: 167)

cgtcgttctgcgaactacgaacctaacagctgggactatgattacctgctgtcctccga
cacggacgagtcctcgaagtatacaaagacaaagcgaaaaagctggaagccgaagttc
gtcgcgagattaataacgaaaaagcagaatttctgaccctgctggaactgattgacaac
gtccagcgcctgggcctgggttaccgtttcgagtctgatccgtcgtgcgctggatcg
cttcgtttcctccggcggcttcgatgcggtaaccaagacttcctgcacgcgacggcac
tgtctttcgtctgctgctcaacacggttttgaggtttctcaggaagcgttcagcggc
ttcaaagacccaaacggcaacttcctggagaacctgaaggaagatatcaaagctatcct
gagcctgtacgagggcagcttcctggctctggaaggcgaaaacatcctggacgaggcga
aggttttcgcaatctctcatctgaaagaactgtctgaagaaaagatcggtaaagatctg
gcagaacaggtgaacctgcaactggaactgccactgcatcgccgtactcagcgtctgga
agcagtactgtctatcgaggcctaccgtaaaaaggaggacgcggatcaggttctgctgg
agctggcaattctggattacaacatgatccagtctgtataccagcgtgatctgctgaa
acgtcccgttgggtggcgtcgtgtgggtctggcgacccaaactgcactttgctcgtgaccg
cctgattgagagcttctactgggcgtgggtgtagcattcgaaccgcaatactccgact
ggcgtaaactccgtcgcaaaaatgttttctttcgttaaccattatcgacgatatctacgat
gtatacggcacccctggacgaactggagctgttactaacgcagttgagcgttgggacgt
aaacgccatcgacgatctgccggattacatgaaactgtgctttctggctctgtataaca
ctattaacgaaatcgctacgacaacctgaaagaaaaaggtgagaacatcctgccgat
ctgaccaaagcctgggctgacctgtgcaacgctttcctgcaagaagccaagtggctgta
caacaaatctactccgacctttgacgaatacttcggcaacgcatggaatcctcttctg
gcccgtgcaactgggtgttcgcttacttcgctgtcgtgcagaacattaaaaaggaagag
atcgaaaacctgcaaaaataccatgacatcatctcgtccttcccatatcttccgtct
gtgcaatgacctggctagcgcgtctgcggaaattgcgcgtggtgaaaccgcaaatagcg
tttcttgttacatgcgcactaaaggtatctccgaagaactggctaccgaaagcgtgatg
aatctgatcgatgaaacctggaaaaagatgaacaaggaaaaactgggtggtagcctgtt
cgcgaaaccgttcgtggaaaccgcgatcaacctggcacgtcaatctcactgcacttatc
ataacggcgacgcgcatacctctccggatgagctgaccgcgaaacgcgttctgtctgta
atcactgaaccgattctgccgtttgaacgctaaaaggcgagctcaacgatccggctgc
taacaaagcccgaaggaagctgagttggctgctgccaccgctgagcaataactagcat
aacccttggggcctctaaacgggtcttgaggagtttttgctgaaaggaggaactata
tccggatatcccgaagaggcccgagcaccggcataaccaagcctatgcctacagca
tccaggggtgacgggtgccgaggatgacgatgagcgcattgttagatttcatacacggtgc
ctgactgcgttagcaatttaactgtgataaaactaccgcattaaagcttatcgatgataa
gctgtcaaacatgagaattaattcttgaagacgaaagggcctcgtgatacgcctatttt
tataggttaatgtcatgataataatggtttcttagacgtcaggtggcacttttcgggga
aatgtgcgcggaacccctatttgtttattttctaaatacattcaaataatgtatccgct
catgagacaataaacctgataaatgcttcaataatattgaaaaaggaagagtatgattg
aacaagatggattgcacgcaggttctccggccgcttgggtggagaggctattcggctat
gactgggcacaaactgacaatcggctgctctgatgccgccgtgttccggctgtcagcgca
ggggcgcccggttctttttgtcaagaccgacctgtccggtgccctgaatgaactgcagg
acgaggcagcgcggctatcgtggctggccacgacggcggttccttgcgcagctgtgctc
gacgttgtcactgaagcgggaagggaactggctgctattgggcgaagtgccggggcagga
tctcctgtcatctcaccttgctcctgccgagaaagtatccatcatggctgatgcaatgc

Figure 90B

Nucleotide Sequence of pDu32

ggcggctgcatacgccttgatccggctacctgccattcgaccaccaagcgaaacatcgc
atcgagcgggcacgtactcggatggaagccggctcttgatcgatcaggatgatctggacga
agagcatcaggggctcgccgagccgaactgttcgccaggctcaaggcgcgcacgccc
acggcgaggatctcgtcgtgacacatggcgatgcctgcttgccgaatatcatggtggaa
aatggccgctttctggattcatcgactgtggccggctgggtgtggcggaccgctatca
ggacatagcgttggctacccgtgatattgctgaagagcttggcggcgaatgggctgacc
gcttcctcgtgctttacggatcgccgctcccgattcgacgcgcacgccttctatcgc
cttcttgacgagttcttctgagcgggactctggggttcgaaatgaccgaccaagcgacg
cctaactgtcagaccaagtttactcatatatactttagattgatttaaaacttcatttt
taatttaaaaggatctaggtgaagatcctttttgataatctcatgaccaaactccctta
acgtgagttttcgttcactgagcgtcagaccccgtagaaaagatcaaaggatcttctt
gagatccttttttctgcgcgtaatctgctgcttgcaaacaaaaaaccaccgctacca
gcggtggtttgtttgccggatcaagagctaccaactcttttcgaaggtaactggctt
cagcagagcgcagataccaaatactgtccttctagtgtagccgtagttaggccaccact
tcaagaactctgtagcaccgcctacatacctcgctctgctaactcctgttaccagtggct
gctgccagtggcgataagtcgtgtcttaccgggttggaactcaagacgatatgttaccgga
taaggcgcagcggctcgggctgaacggggggttcgtgcacacagcccagcttgagcgaa
cgacctacaccgaactgagatacctacagcgtgagctatgagaaagcgcacgcttccc
gaaggagaaaaggcggacaggtatccggtgaagcggcagggtcggaacaggagagcgac
gaggagacttccaggggaaacgcctggtatctttatagtccctgtcgggtttcgccacc
tctgacttgagcgtcgattttttgtgatgctcgtcagggggcgagcctatggaaaaac
gccagcaacgcggcctttttacgggttcctggccttttgctggccttttgctcacatgtt
ctttcctgcgttatccctgattctgtggataaccgtattaccgcctttgagtgaagctg
ataccgctcgccgcagccgaacgaccgagcgcagcagtcagtgaagcaggaaagcggaa
gagcgcctgatgcggatattttctccttacgcacatctgtgcggatattcacaccgcaatgg
tgcaactctcagtaacaatctgctctgatgccgcatagttaagccagtatacactccgcta
tcgctacgtgactgggtcatggctgcgccccgacaccgcgaacaccgcgtgacgcgcc
ctgacgggcttgtctgctcccggcatccgcttacagacaagctgtgaccgtctccggga
gctgcatgtgtcagaggttttaccgctcatcaccgaaacgcgcgaggcagctgcggtaa
agctcatcagcgtggctcgtgaagcgattcacagatgtctgcctgttcatccgcgtccag
ctcgttgagttttctccagaagcgttaatgtctggcttctgataaagcgggccaatgttaa
gggcggttttttctggttggctcactgatgcctccgtgaaggggatttctgttcatg
ggggtaatgataccgatgaaacgagagaggatgctcacgatacgggttactgatgatga
acatgccgggttactggaacgttgtgagggttaacaactggcggatggatgcgcggg
accagagaaaaatcactcagggccaatgccagcgttcgtaataacagatgtaggtgtt
ccacagggtagccagcagcatcctgcgatgcagatccggaacataatggtgcagggcgc
tgacttccgcgtttccagactttacgaaacacgaaaccgaagaccattcatgttgttg
ctcaggtcgacagcgttttgacgagcagtcgcttcacgttcgctcgctatcggtgat
tcattctgctaaccagtaaggcaaccccgccagcctagccgggtcctcaacgacaggag
cacgatcatgcgcacccgtggccaggaccaacgctgcccgagatgcgcgcgtgcggc
tgctggagatggcggacgcgatggatatgttctgccaagggttggtttgcgcattcaca
gttctccgcaagaattgattggctccaattcttgagtggtgaatccgttagcagagtg
ccgcgggcttccattcaggtcgaggtggcccggtccatgcaccgcgacgcaacgcggg

Figure 90C

Nucleotide Sequence of pDu32

gaggcagacaaggtatagggcggcgcctacaatccatgccaaaccggttccatgtgctcg
ccgagggcggcataaatcgccgtgacgatcagcgggtccaatgatcgaagttaggctggta
agagccgcgagcgatccttgaagctgtccctgatggtcgtcatctacctgcctggacag
catggcctgcaacgcgggcatcccgatgccgcggaagcgagaagaatcataatgggga
aggccatccagcctcgctcggaacgccagcaagacgtagcccagcgcgtcgcccgcc
atgccggcgataatggcctgcttctcgccgaaacgtttggtggcgggaccagtgcgaa
ggcttgagcgagggcggtgcaagattccgaataccgcaagcgacaggccgatcatcgtcg
cgctccagcgaaagcggtcctcgccgaaaatgaccagagcgctgccggcacctgtcct
acgagttgcatgataaagaagacagtcataagtgcggcgacgatagtcagccccgcgc
ccaccggaaggagctgactgggttgaaggctctcaagggcacggctcgagatcccggtg
cctaagtagtgaactaacattaatgagcttgcgttgcgtcactgcccgtttccagtcg
ggaaacctgtcgtgccagctgcattaatgaatcgcccaacgcgcggggagaggcggttt
gcgtattgggcgccagggtgggtttttcttttaccagtgagacgggcaacagctgattg
cccttcaccgcctggccctgagagagttgcagcaagcggtccacgctggtttgccccag
caggcgaaaatcctgtttgatgggtgggttaacggcgggatataacatgagctgtctcgg
tatcgtcgtatcccactaccgagatccgcaccaacgcgcagcccggactcggtaatg
gcgcgcatttgcgccagcgccatctgatcgttggcaaccagcatcgagtggaacgat
gccctcattcagcatttgcattggttgggtgaaaaccggacatggcactccagtcgcctt
cccggttccgctatcggtgaatttgattgcgagtgagatatttatgccagccagccaga
cgcagacgcgcgcgagacagaacttaatgggcccgcctaacagcgcgatttgctggtgacc
caatgcgaccagatgctccacgcccagtcgcgtaccgtcttcatgggagaaaataatac
tgttgatgggtgtctgggtcagagacatcaagaaataacgcccgaacattagtgcaggca
gcttccacagcaatggcatcctggtcatccagcggatagttaatgatcagcccactgac
gcgttgcgcgagaagattgtgcaccgcgcgtttacaggcttcgacgcgcgttcgttcta
ccatcgacaccaccagctggcaccagttgatcggcgcgagatttaatcgccgcgaca
atttgcgacggcgcggtgcagggccagactggaggtggcaacgccaatcagcaacgactg
tttgcccgcagttgttgcacgcgggttggaatgtaattcagctccgccatcgccg
cttccactttttcccgcttttcgcagaaacgtggctggcctggttcaccacgcgggaa
acggtctgataagagacaccggcatactctgcgacatcgtataacgttactggtttcac
attcaccaccctgaattgactctcttccggcgctatcatgccataccgcgaaaggttt
tgcgccattcgatgggtgtccgggatctcgacgctctcccttatgcgactcctgcattag
gaagcagcccagtagtaggttgaggccgttgagcaccgcgcgcgcaaggaatggtgcat
gcaaggagatggcgcccaacagtcggcgccacggggcctgccaccataccacgccc
aaacaagcgctcatgagcccgaagtggcgagcccgatcttccccatcggtgatgtcggc
gatataggcgccagcaaccgcacctgtggcgccggtgatgccggccacgatgcgtccgg
cgtagaggatcgagatctcgatcccgcgaaattaatacgactcactataggggaattgt
gagcggataacaattcccctctagaaataattttgtttaactttaagaaggagatatc
atatgccccgttctcatcatcatcatcatcatggtatggctagcatgactggtggacag
caaatgggtcgggatctgtacgacgatgacgataaggatcatcccttcaccatg

Figure 91




Figure 92




Figure 93




Figure 94




Figure 95




Figure 96




Figure 97

