

(86) Date de dépôt PCT/PCT Filing Date: 2012/03/30

(87) Date publication PCT/PCT Publication Date: 2012/10/04

(45) Date de délivrance/Issue Date: 2019/05/21

(85) Entrée phase nationale/National Entry: 2013/09/27

(86) N° demande PCT/PCT Application No.: US 2012/031553

(87) N° publication PCT/PCT Publication No.: 2012/135682

(30) **Priorité/Priority:** 2011/03/31 (US61/470,275)

(51) Cl.Int./Int.Cl. *H01J 49/04* (2006.01)

(72) Inventeurs/Inventors:

DEY, SUBHAKAR N., US;
MINKOFF, MARJORIE S., US;
PURKAYASTHA, SUBHASISH, US;
PILLAI, SASI K., US;
WILLIAMSON, BRIAN L., US


(73) Propriétaire/Owner:

DH TECHNOLOGIES DEVELOPMENT PTE. LTD., SG

(74) **Agent:** PERRY + CURRIER

(54) Titre : COMPOSITION, PROCEDE ET NECESSAIRE POUR ETALONNAGE DE SPECTROMETRE DE MASSE

(54) Title: COMPOSITION, METHOD, AND KIT FOR CALIBRATING A MASS SPECTROMETER



(57) **Abrégé/Abstract:**

A composition, method, and kit for calibrating a mass spectrometer including a predetermined concentration of a calibrant having a mixture of amino acid polyethylene glycol compounds and a solvent for dissolving the calibrant. The calibrant can be used in either positive or negative ionization mode, and it can be used in calibrating an atmospheric pressure chemical ionization or an electrospray mass spectrometer. The calibrant can include a mixture to enable calibration across a broad range of masses.

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)
(19) World Intellectual Property Organization
International Bureau

(43) International Publication Date
4 October 2012 (04.10.2012)



(10) International Publication Number
WO 2012/135682 A2

- (51) International Patent Classification:
H01J 49/04 (2006.01) *G01N 27/62* (2006.01)
- (21) International Application Number:
PCT/US2012/031553
- (22) International Filing Date:
30 March 2012 (30.03.2012)
- (25) Filing Language: English
- (26) Publication Language: English
- (30) Priority Data:
61/470,275 31 March 2011 (31.03.2011) US
- (71) Applicant (for all designated States except US): **DH TECHNOLOGIES DEVELOPMENT PTE. LTD.** [SG/SG]; C/o KW Corporate Advisory Pte. Ltd, 80 Raffles Place #25-01, UOB Plaza, Singapore 0486624 (SG).
- (71) Applicants (for US only): **DEY, Subhakar, N.** [IN/US]; 9 Sewell Street, North Billerica, MA 01862 (US). **MINKOFF, Marjorie, S.** [US/US]; 127 Freeman Street, Apt. 3a, Brookline, MA 02446 (US). **PILLAI, Sasi, K.** [US/US]; 1 Nicole Lane, Littleton, MA 01460 (US). **PURKAYASTHA, Subhasish** [IN/US]; 19 Squirrel Hill

Road, Acton, MA 01720 (US). **WILLIAMSON, Brian, L.** [GB/US]; 118 Olive Street, Ashland, MA 01721 (US).

- (81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IS, JP, KE, KG, KM, KN, KP, KR, KZ, LA, LC, LK, LR, LS, LT, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.
- (84) Designated States (unless otherwise indicated, for every kind of regional protection available): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

[Continued on next page]

(54) Title: COMPOSITION, METHOD, AND KIT FOR CALIBRATING A MASS SPECTROMETER

(57) Abstract: A composition, method, and kit for calibrating a mass spectrometer including a predetermined concentration of a calibrant having a mixture of amino acid polyethylene glycol compounds and a solvent for dissolving the calibrant. The calibrant can be used in either positive or negative ionization mode, and it can be used in calibrating an atmospheric pressure chemical ionization or an electrospray mass spectrometer. The calibrant can include a mixture to enable calibration across a broad range of masses.




Figure 1

WO 2012/135682 A2

WO 2012/135682 A2



Published:

- *without international search report and to be republished
upon receipt of that report (Rule 48.2(g))*

**TITLE: COMPOSITION, METHOD, AND KIT FOR CALIBRATING
A MASS SPECTROMETER**

5 RELATED APPLICATION

[0001] This application claims priority to US provisional application no. 61/470,275 filed March 31, 2011.

10 FIELD

[0002] The applicants' teachings relate to a composition, method, and kit for calibrating a mass spectrometer.

INTRODUCTION

15 [0003] To achieve accuracy and reliability, mass spectrometers need to be calibrated. Although there are a variety of calibrants used to calibrate mass spectrometers, many of them do not span a large mass range, do not work in both positive and negative ion modes, and are not ionizable in both electrospray (ESI) and atmospheric pressure ionization (APCI) modes.
20 Typically, multiple sets of compounds are used that work either in positive or negative ion mode and either in ESI or APCI mode.

SUMMARY

[0004] In accordance with an aspect of the applicants' teachings, a calibration composition for use in mass spectrometry is provided. The
25 system comprises a predetermined concentration of a calibrant comprising a mixture of amino acid polyethylene glycol compounds, also known as discrete polyethylene glycol (dPEG[®]) compounds, and a solvent for dissolving the calibrant. In various embodiments, the calibrant can be used in either positive or negative ionization mode. In various embodiments, the calibrant can be
30 used in calibrating an atmospheric pressure chemical ionization (APCI) or an

- 2 -

electrospray mass spectrometer. In various embodiments, the calibrant composition comprises a mixture to enable calibration across a range of approximately 145 to 3500 Da. In various aspects, in calibrating an APCI mass spectrometer in positive ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In various aspects, in calibrating an APCI mass spectrometer in negative ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In various embodiments, in calibrating an electrospray mass spectrometer in positive ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In various embodiments, in calibrating an electrospray mass spectrometer in negative ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used.

[0005] In various aspects, a method for calibrating a mass spectrometer is provided. The method comprises obtaining a mass spectrum

- 3 -

of a calibrant composition containing a plurality of known compounds comprising a mixture of amino acid polyethylene glycol compounds, also known as discrete polyethylene glycol (dPEG[®]) compounds, and a solvent for dissolving the calibrant, determining the differences between the expected
5 mass peaks for the known compounds and the corresponding actual mass peaks obtained, and adjusting the mass spectrometer based on the differences between the expected and actual mass peaks. In various embodiments, the calibrant can be used in either positive or negative ionization mode. In various embodiments, the calibrant can be used in
10 calibrating an atmospheric pressure chemical ionization (APCI) or an electrospray mass spectrometer. In various embodiments, the calibrant mixture is selected to enable calibration across a range of approximately 145 to 3500 Da. In various aspects, in calibrating an APCI mass spectrometer in positive ionization mode, the polyethylene glycol compounds comprise amino-
15 dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In various
20 aspects, in calibrating an APCI mass spectrometer in negative ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone. In various aspects, the solvent
25 comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In various embodiments, in calibrating an electrospray mass spectrometer in positive ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the calibrant further
30 comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In

- 4 -

various embodiments, in calibrating an electrospray mass spectrometer in negative ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various embodiments, the calibrant
5 further comprises 7-aminoheptanoic acid and sulfinpyrazone. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used.

[0006] In various aspects, a kit for calibrating a mass spectrometer is provided comprising a predetermined concentration of a
10 calibrant comprising a mixture of amino acid polyethylene glycol compounds, also known as discrete polyethylene glycol (dPEG) compounds, and a solvent for dissolving the calibrant. In various embodiments, the calibrant comprises a mixture to enable calibration across a range of approximately 145 to 3500 Da. In various aspects, in calibrating an atmospheric pressure chemical
15 ionization (APCI) mass spectrometer in positive ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521. In
20 various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In various aspects, in calibrating an atmospheric pressure chemical ionization (APCI) mass spectrometer in negative ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG
25 12-acid, and amino-dPEG 16-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In various embodiments, in calibrating an electrospray mass spectrometer in positive ionization mode, the polyethylene
30 glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the

- 5 -

calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521. In various aspects, the solvent comprises a mixture of acetonitrile and water, but any suitable solvent can be used. In various embodiments, in calibrating an electrospray mass spectrometer in negative ionization mode, the polyethylene glycol compounds
5 comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone. In various aspects, the solvent comprises a mixture of acetonitrile and water, but
10 any suitable solvent can be used.

[0007] These and other features of the applicants' teachings are set forth herein.

BRIEF DESCRIPTION OF THE DRAWINGS

[0008] The skilled person in the art will understand that the
15 drawings, described below, are for illustration purposes only. The drawings are not intended to limit the scope of the applicants' teachings in anyway.

[0009] Figure 1 shows the molecular structure of discrete polyethylene compounds in accordance with various embodiments of the applicants' teachings.

20 [0010] Figure 2 shows a table of a calibration composition for an APCI ion source in positive ionization mode in accordance with various embodiments of the applicants' teachings.

[0011] Figure 3 shows mass spectra of a calibrant composition obtained using a mass spectrometer with an APCI ion source in positive
25 ionization mode in accordance with various embodiments of the applicants' teachings.

- 6 -

[0012] Figure 4 shows an MS/MS mass spectra of clomipramine obtained with an APCI source in positive ionization mode in accordance with various embodiments of the applicants' teachings.

[0013] Figure 5 shows an MS/MS mass spectra of reserpine
5 obtained with an APCI source in positive ionization mode in accordance with various embodiments of the applicants' teachings.

[0014] Figure 6 shows an MS/MS mass spectra of verapamil obtained with an APCI source in positive ionization mode in accordance with various embodiments of the applicants' teachings.

10 [0015] Figure 7 shows mass spectra of various compounds after calibrating with a calibrant composition using an APCI ion source in positive ionization mode in accordance with various embodiments of the applicants' teachings.

[0016] Figure 8 shows a table of a calibration composition for an
15 APCI ion source in negative ionization mode in accordance with various embodiments of the applicants' teachings.

[0017] Figure 9 shows mass spectra of a calibrant composition obtained using a mass spectrometer with an APCI ion source in negative ionization mode in accordance with various embodiments of the applicants'
20 teachings.

[0018] Figure 10 shows an MS/MS mass spectra of sulfinpyrazone obtained with an APCI source in negative ionization mode in accordance with various embodiments of the applicants' teachings.

[0019] Figure 11 shows an MS/MS mass spectra of
25 bromocriptine obtained with an APCI source in negative ionization mode in accordance with various embodiments of the applicants' teachings.

- 7 -

[0020] Figure 12 shows mass spectra of Bromocriptine after calibrating with a calibrant composition using an APCI ion source in negative ionization mode in accordance with various embodiments of the applicants' teachings.

5 [0021] Figure 13 shows a table of a calibration composition for an ESI ion source in positive ionization mode in accordance with various embodiments of the applicants' teachings.

[0022] Figure 14 shows mass spectra of a calibrant composition obtained using a mass spectrometer with an ESI ion source in positive
10 ionization mode in accordance with various embodiments of the applicants' teachings.

[0023] Figure 15 shows mass spectra of a calibrant composition obtained using a mass spectrometer with an ESI ion source in positive
15 ionization mode in accordance with various embodiments of the applicants' teachings.

[0024] Figure 16 shows a table of a calibration composition for an APCI ion source in negative ionization mode in accordance with various embodiments of the applicants' teachings.

[0025] Figure 17 shows mass spectra of a calibrant composition
20 obtained using a mass spectrometer with an ESI ion source in negative ionization mode in accordance with various embodiments of the applicants' teachings.

[0026] Figure 18 shows mass spectra of a calibrant composition obtained using a mass spectrometer with an ESI ion source in negative
25 ionization mode in accordance with various embodiments of the applicants' teachings.

- 8 -

[0027] Figure 19 illustrates the stability of a calibrant composition in positive ionization mode in accordance with various embodiments of the applicants' teachings.

[0028] Figure 20 illustrates the stability of a calibrant composition
5 in negative ionization mode in accordance with various embodiments of the applicants' teachings.

DESCRIPTION OF VARIOUS EMBODIMENTS

[0029] It should be understood that the phrase "a" or "an" used in conjunction with the applicants' teachings with reference to various elements
10 encompasses "one or more" or "at least one" unless the context clearly indicates otherwise. In various embodiments, a calibration composition comprises a predetermined concentration of a calibrant comprising a mixture of amino acid polyethylene glycol compounds or discrete polyethylene glycol (dPEG[®]) compounds in a solvent that can dissolve the calibrant and any other
15 components that can be present in the composition. In various embodiments, the calibrant composition can be used in positive or negative ionization mode. In various aspects, certain calibrant compositions can be better suited for use in positive ionization mode while others can be better for use in negative ionization mode. In various embodiments, the calibrant composition can be
20 used in calibrating a mass spectrometer with an atmospheric pressure chemical ionization (APCI) mass spectrometer or a mass spectrometer with an electrospray (ESI) ion source. In various embodiments, a calibrant composition can be selected to enable calibration across a broad mass range and can provide reference mass peaks in both positive and negative
25 ionization modes. In various aspects, substantially the same components can be used for both APCI and electrospray mass spectrometry. To provide calibration across a broad mass range, the calibrant composition can comprise a mixture of different discrete polyethylene glycol compounds in accordance with various embodiments of the applicant's teachings. In various
30 aspects, the calibrant composition can enable calibration across a range of

- 9 -

approximately 145 to 3500 Da. In various aspects, the solvent to dissolve the calibrant can comprise a mixture of acetonitrile and water, but any suitable solvent can be used.

[0030] Reference is made to Figure 1 which shows the molecular structure of discrete polyethylene glycol compounds that can comprise the calibrant composition in accordance with various embodiments of the applicants' teachings. In various aspects, the calibrant composition can include, but is not limited to, the discrete polyethylene glycol compounds shown in Figure 1. In various embodiments, as shown in Figure 1, the calibration composition can comprise a mixture of discrete length polyethylene glycol compounds having an amino group on one end, a number of ethylene oxide units, and a carboxylic acid group on the other end, known also as amino-dPEG_n-acids, wherein n is the number of ethylene oxide units. For example, amino-dPEG₄-acid has an amino group on one end, four ethylene oxide units, and a carboxylic acid group on the other end, as shown in Figure 1. In various embodiments, n can be in the range of 4 to 16. In various aspects, the calibrant composition comprises a mixture of different amino-dPEG-acid compounds that can be selected to span across a broad mass range. In various embodiments, approximately four to five amino-dPEG-acid compounds can be used in the calibration composition. For example, the following amino-dPEG-acid compounds in the calibration composition can include, but are not limited to, amino-dPEG₄-acid, amino-dPEG₆-acid, amino-dPEG₈-acid, amino-dPEG₁₂-acid, and amino-dPEG₁₆-acid. The amino-dPEG-acid compounds can be obtained commercially from QuantaBiodesign. In various aspects, the amino-dPEG_n-acid compounds in the calibration composition can be produced by mixing an assortment of amino-dPEG_n and acid-dPEG_n of discrete masses. The calibration composition can further include, but is not limited to, 7-aminoheptanoic acid (which can be commercially obtained from Sigma-Aldrich), clomipramine (which can be commercially obtained from Sigma-Aldrich), reserpine (which can be commercially obtained from Sigma-Aldrich), phosphazene 921 (which can be

- 10 -

commercially obtained from Apollo Scientific Ltd.), phosphazene, 1521 (which can be commercially obtained from Apollo Scientific Ltd.), and sulfinpyrazone (which can be commercially obtained from Sigma-Aldrich).

[0031] In various embodiments, in calibrating an APCI mass spectrometer in positive ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.

10 [0032] In various embodiments, a calibration composition for use in APCI mass spectrometry calibration in positive ionization mode was prepared comprising the components in Figure 2. Figure 3 shows the mass spectra of the calibration composition obtained in APCI positive ionization mode from two instruments according to various embodiments of the
15 applicant's teachings. In various aspects, Figure 4 shows an MS/MS spectrum obtained for clomipramine in positive ionization mode following an autotune calibration with the APCI calibration composition prepared in Figure 2. The individual mass peaks found closely agreed with the target mass peaks as shown in the table in Figure 4. Figure 5 shows an MS/MS spectrum
20 obtained for reserpine in positive ionization mode following an autotune calibration with the APCI calibration composition prepared in Figure 2. The individual mass peaks found closely agreed with the target mass peaks as shown in the table in Figure 5. Figure 6 shows an MS/MS spectrum obtained for verapamil in positive ionization mode following an autotune calibration with
25 the APCI calibration composition prepared in Figure 2. The individual mass peaks found closely agreed with the target mass peaks as shown in the table in Figure 6. According to various embodiments of the applicant's teachings, Figure 7 shows the mass accuracy obtained for various compounds in positive ionization mode following an autotune calibration with the APCI
30 calibration composition prepared in Figure 2.

- 11 -

[0033] In various aspects, in calibrating an APCI mass spectrometer in negative ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various embodiments, the
5 calibrant composition further comprises 7-aminoheptanoic acid and sulfinpyrazone.

[0034] In various embodiments, a calibration composition for use in APCI mass spectrometry calibration in negative ionization mode was prepared comprising the components in Figure 8. Figure 9 shows the mass
10 spectra of the calibration composition obtained in APCI negative ionization mode from two instruments according to various embodiments of the applicant's teachings. In various aspects, Figure 10 shows an MS/MS spectrum obtained for sulfinpyrazone in negative ionization mode following an autotune calibration with the APCI calibration composition prepared in Figure
15 8. The individual mass peaks found closely agreed with the target mass peaks as shown in the table in Figure 10. Figure 11 shows an MS/MS spectrum obtained for bromocriptine in negative ionization mode following an autotune calibration with the APCI calibration composition prepared in Figure
20 8. The individual mass peaks found closely agreed with the target mass peaks as shown in Figure 11. According to various embodiments of the applicant's teachings, Figure 12 shows the mass accuracy obtained for bromocriptine in negative ionization mode following an autotune calibration with the APCI calibration composition prepared in Figure 8.

[0035] In various embodiments, in calibrating an electrospray
25 mass spectrometer in positive ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid clomipramine, reserpine, phosphazene 921, and phosphazene 1521.

- 12 -

[0036] In various embodiments, a calibration composition for use in electrospray mass spectrometry calibration in positive ionization mode was prepared comprising the components in Figure 13. Figure 14 shows the mass spectrum of the calibration composition obtained in electrospray positive ionization mode according to various embodiments of the applicant's teachings. In various aspects, Figure 15 shows the mass spectrum obtained from a 100 ng/mL mixture of calibration composition in electrospray positive ionization mode.

[0037] In various embodiments, in calibrating an electrospray mass spectrometer in negative ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various embodiments, the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone.

[0038] In various embodiments, a calibration composition for use in electrospray mass spectrometry calibration in negative ionization mode was prepared comprising the components in Figure 16. Figure 17 shows the mass spectrum of the calibration composition obtained in electrospray negative ionization mode according to various embodiments of the applicant's teachings. Figure 18 shows a spectrum obtained in negative ionization mode following an autotune calibration with the ESI calibration composition prepared in Figure 16. The individual mass peaks found closely agreed with the target mass peaks as shown in the table in Figure 18.

[0039] According to various embodiments of the applicant's teachings, Figure 19 shows the stability of the positive ion calibration solution during a span of five weeks at a temperature of 25 degrees Celsius. Figure 20 shows the stability of the negative ion calibration solution during a span of five weeks at a temperature of 25 degrees Celsius according to various embodiments of the applicant's teachings.

- 13 -

[0040] In various embodiments, solvents can be used to dissolve the calibrant and can include, but are not limited to, acetonitrile, methanol, ethanol, propanol, isopropanol, or water. In various aspects, the solvent can comprise a mixture of acetonitrile and water.

5 **[0041]** In various aspects, the calibration composition can be used as an internal standard to optimize and tune parameters of the mass spectrometer to ensure that accurate mass spectra are obtained. In use, a mass spectrum of a calibrant composition can be obtained. The calibrant composition can include a plurality of known compounds comprising a mixture
10 of amino acid polyethylene glycol compounds and a solvent for dissolving the calibrant. The differences between the expected mass peaks for the known compounds and the corresponding mass peaks obtained can be determined and the mass spectrometer can be adjusted based on the differences between the expected and actual mass peaks. In various aspects, the
15 calibrant can be used in either positive or negative ionization mode, and the calibrant can be used in calibrating an APCI or an electrospray mass spectrometer.

[0042] In various embodiments, in calibrating an APCI mass spectrometer in positive ionization mode, the polyethylene glycol compounds
20 can comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the calibrant can further comprise 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.

[0043] In various embodiments, in calibrating an APCI mass
25 spectrometer in negative ionization mode, the polyethylene glycol compounds can comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various aspects, the calibrant can further comprise 7-aminoheptanoic acid and sulfinpyrazone.

- 14 -

[0044] In various embodiments, in calibrating an electrospray mass spectrometer in positive ionization mode, the polyethylene glycol compounds can comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various aspects, the calibrant can
5 further comprise 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.

[0045] In various embodiments, in calibrating an electrospray mass spectrometer in negative ionization mode, the polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG
10 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various aspects, the calibrant can further comprise 7-aminoheptanoic acid and sulfinpyrazone.

[0046] In various embodiments, solvents can be used to dissolve the calibrant and can include, but are not limited to, acetonitrile, methanol, ethanol, propanol, isopropanol, or water. In various aspects, the solvent can
15 comprise a mixture of acetonitrile and water.

[0047] In various embodiments, the calibrant composition can include a mixture to enable calibration across a broad range of masses. In various aspects, the calibrant composition can enable calibration across a range of approximately 145 to 3500 Da.

20 **[0048]** In various embodiments, a kit for calibrating a mass spectrometer can be provided. The kit can comprise a predetermined concentration of a calibrant comprising a mixture of amino acid polyethylene glycol compounds and a solvent for dissolving the calibrant.

[0049] In various aspects, in calibrating an APCI mass
25 spectrometer in positive ionization mode, the polyethylene glycol compounds in the kit can comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various embodiments, the kit can further comprise 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.

- 15 -

[0050] In various embodiments, in calibrating an APCI mass spectrometer in negative ionization mode, the polyethylene glycol compounds in the kit can comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various aspects, the kit can further comprise 7-aminoheptanoic acid and sulfinpyrazone.

[0051] In various embodiments, in calibrating an electrospray mass spectrometer in positive ionization mode, the polyethylene glycol compounds in the kit can comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid. In various aspects, the kit can further comprise 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.

[0052] In various embodiments, in calibrating an electrospray mass spectrometer in negative ionization mode, the polyethylene glycol compounds in the kit can comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid. In various aspects, the kit can further comprise 7-aminoheptanoic acid and sulfinpyrazone.

[0053] In various aspects, the kit can comprise a solvent to dissolve the calibrant. In various embodiments, the solvent can comprise a mixture of acetonitrile and water. In various embodiments, the kit can contain any suitable solvent.

[0054] In various embodiments, the calibrant composition can include a mixture to enable calibration across a broad range of masses. In various aspects, the kit can comprise a calibrant composition that can enable calibration across a range of approximately 145 to 3500 Da.

[0055] While the applicants' teachings are described in conjunction with various embodiments, it is not intended that the applicants' teachings be limited to such embodiments. On the contrary, the applicants' teachings encompass various alternatives, modifications, and equivalents, as

- 16 -

will be appreciated by those skilled in the art.

CLAIMS

1. A calibration composition for use in mass spectrometry comprising:
 - 5 a predetermined concentration of a calibrant comprising a mixture of discrete length amino acid polyethylene glycol compounds, the discrete length amino acid polyethylene glycol compounds each having an amino group on one end, a number of ethylene oxide units and a carboxylic acid group on the other end, and;
 - a solvent for dissolving the calibrant.
- 10 2. The composition of claim 1 wherein the calibrant is used in either positive or negative ionization mode.
3. The composition of claim 2 wherein the calibrant is used in calibrating an atmospheric pressure chemical ionization or an electrospray mass spectrometer.
- 15 4. The composition of claim 3 wherein in calibrating an atmospheric pressure chemical ionization mass spectrometer in positive ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid.
- 20 5. The composition of claim 3 wherein in calibrating an atmospheric pressure chemical ionization mass spectrometer in negative ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid.
- 25 6. The composition of claim 4 wherein the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.

7. The composition of claim 5 wherein the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone.
8. The composition of claim 3 wherein in calibrating an electrospray mass spectrometer in positive ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid.
9. The composition of claim 3 wherein in calibrating an electrospray mass spectrometer in negative ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid.
10. The composition of claim 8 wherein the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.
11. The composition of claim 9 wherein the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone.
12. The composition of claim 1 wherein the solvent comprises a mixture of acetonitrile and water.
13. The composition of claim 1 wherein the calibrant comprises a mixture to enable calibration across a range of approximately 145 to 3500 Da.
14. A method for calibrating a mass spectrometer comprising:
 - a) obtaining a mass spectrum of a calibrant composition containing a plurality of known compounds comprising a mixture of discrete length amino acid polyethylene glycol compounds and a solvent for dissolving the calibrant, each of the discrete length amino acid polyethylene glycol compounds having an amino group on one end, a number of ethylene oxide units, and a carboxylic acid group on the other end;
 - b) determining the differences between the expected mass peaks for the known compounds and the corresponding actual mass peaks obtained.

c) adjusting the mass spectrometer based on the differences between the expected and actual mass peaks.

15. The method of claim 14 wherein the calibrant is used in either positive or negative ionization mode.

5 16. The method of claim 15 wherein the calibrant is used in calibrating an atmospheric pressure chemical ionization or an electrospray mass spectrometer.

10 17. The method of claim 16 wherein in calibrating an atmospheric pressure chemical ionization mass spectrometer in positive ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid.

15 18. The method of claim 16 wherein in calibrating an atmospheric pressure chemical ionization mass spectrometer in negative ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid.

20 19. The method of claim 17 wherein the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.

20. The method of claim 18 wherein the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone.

25 21. The method of claim 16 wherein in calibrating an electrospray mass spectrometer in positive ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid.

22. The method of claim 16 wherein in calibrating an electrospray mass spectrometer in negative ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-

dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid.

23. The method of claim 21 wherein the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and
5 phosphazene 1521.

24. The method of claim 22 wherein the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone.

25. The method of claim 14 wherein the solvent comprises a mixture of acetonitrile and water.

10 26. The method of claim 14 wherein the calibrant mixture is selected to enable calibration across a range of approximately 145 to 3500 Da.

27. A kit for calibrating a mass spectrometer comprising:

15 a predetermined concentration of a calibrant comprising a mixture of discrete length amino acid polyethylene glycol compounds, each of the discrete length amino acid polyethylene glycol compounds having an amino acid group on one end, a number of ethylene oxide units, and a carboxylic acid group on the other end, and;

a solvent for dissolving the calibrant.

28. The kit of claim 27 wherein in calibrating an atmospheric pressure
20 chemical ionization mass spectrometer in positive ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid.

29. The kit of claim 27 wherein in calibrating an atmospheric pressure
25 chemical ionization mass spectrometer in negative ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid.

30. The kit of claim 28 wherein the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.
- 5 31. The kit of claim 29 wherein the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone.
32. The kit of claim 27 wherein in calibrating an electrospray mass spectrometer in positive ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, and amino-dPEG 12-acid.
- 10 33. The kit of claim 27 wherein in calibrating an electrospray mass spectrometer in negative ionization mode, the discrete length amino acid polyethylene glycol compounds comprise amino-dPEG 4-acid, amino-dPEG 6-acid, amino-dPEG 8-acid, amino-dPEG 12-acid, and amino-dPEG 16-acid.
- 15 34. The kit of claim 32 wherein the calibrant further comprises 7-aminoheptanoic acid, clomipramine, reserpine, phosphazene 921, and phosphazene 1521.
35. The kit of claim 33 wherein the calibrant further comprises 7-aminoheptanoic acid and sulfinpyrazone.
- 20 36. The kit of claim 27 wherein the solvent comprises a mixture of acetonitrile and water.
37. The kit of claim 27 wherein the calibrant comprises a mixture to enable calibration across a range of approximately 145 to 3500 Da.




Figure 1

APCI Positive Calibration Solution									
Name	Formula (M)	CAS	Sigma	QuantaBiodesign	Apollo	FW	molarity		
7-aminoheptanoic acid	C7H15NO2	929-17-9	284637			145.2	6.89E-06		
amino-dPEG 4-acid	C11H23NO6			10244		265.3	5.65E-07		
clomipramine	C19H23ClN2.HCl	17321-77-6	C7291			351.3 1	4.98E-07		
amino-dPEG 6-acid	C15H31NO8			10067		353.4 1	2.12E-07		
amino-dPEG 8-acid	C18H39NO10	756526-04-2		10277		441.5 1	1.93E-07		
Reserpine	C33H40N2O9	50-55-5	R0875			608.6 8	6.24E-07		
amino-dPEG 12-acid	C27H55NO14	756526-07-5		10287		617.7 2	4.61E-07		
phosphazene 921	C18H18F24N3O6P3	58943-98-9			PC087 4	921.2 3	3.05E-07		
phosphazene 1521	C30H18F48N3O6P3	16059-16-8			PC023 3	1521. 33	3.29E-07		

Figure 2

Positive Ion Spectra (APCI)





Figure 3

Positive MS/MS - clomipramine

■ +TOF Product (315.1): Exp 2, 3.4576 to 3.4880 min from Sample 2 (p9 ext-002) of P9 ext cal 1-10, wiff different calibrations (DuoSpray (I)) Max. 595.3 cps.



Peak List for +TOF Product (315.1): Exp 2, 3.4576 to 3.4880 min from Sample 2 (p9 ext-002) of P9 ext cal 1-10, wiff different calibrations (DuoSpray (I))				
	Target Mass (Da)	Found At (Da)	Intensity (cps)	Width (Da)
1	58.06510	58.06514	229.75000	3.46332e-3
2	86.09640	86.09669	595.25000	4.51873e-3
3	242.07310	242.07333	163.25000	6.59612e-3
4	270.10440	270.10442	74.00000	6.52390e-3
5	315.16225	315.16281	301.25000	0.01107

Figure 4

Positive MS/MS - reserpine




Figure 5

Positive MS/MS - verapamil




Figure 6

Mass Accuracy after Calibrating - positive





Figure 7

APCI Negative Calibration Solution									
Name	Formula (M)	CAS	Sigma	QuantaBiodesign	FW	molarity			
7-aminoheptanoic acid	C7H15NO2	929-17-9	284637		145.2	6.89E-06			
amino-dPEG 4-acid	C11H23NO6			10244	265.3	9.42E-07			
amino-dPEG 6-acid	C15H31NO8			10067	353.41	7.07E-07			
sulfinpyrazone	C23H20N2O3S	57-96-5	S9509		404.48	8.65E-07			
amino-dPEG 8-acid	C18H39NO10	756526-04-2		10277	441.51	2.83E-07			
amino-dPEG 12-acid	C27H55NO14	756526-07-5		10287	617.72	4.61E-07			
amino-dPEG 16-acid	C35H71NO18			10297	793.93	7.87E-07			

Figure 8

Negative Ion Spectra (APCI)



Negative MS/MS Calibration - sulfinpyrazone




Figure 10

Negative MS/MS - bromocriptine




Figure 11

**T60
Bromocriptine incubation**



Figure 12

ESI Positive Calibration Solution										
Name	Formula (M)	CAS	Sigma	QuantaBiodesign	Apollo	FW	molarity			
7-aminoheptanoic acid	C7H15NO2	929-17-9	284637			145.2	3.44E-06			
amino-dPEG 4-acid	C11H23NO6			10244		265.3	2.83E-07			
clomipramine	C19H23ClN2.HCl	17321-77-6	C7291			351.31	2.49E-07			
amino-dPEG 6-acid	C15H31NO8			10067		353.41	1.06E-07			
amino-dPEG 8-acid	C18H39NO10	756526-04-2		10277		441.51	9.63E-08			
Reserpine	C33H40N2O9	50-55-5	R0875			608.68	1.25E-06			
amino-dPEG 12-acid	C27H55NO14	756526-07-5		10287		617.72	1.15E-08			
phosphazene 921	C18H18F24N3O6P3	58943-98-9			PC0874	921.23	3.05E-06			
phosphazene 1521	C30H18F48N3O6P3	16059-16-8			PC0233	1521.33	3.29E-06			

Figure 13

ESI Positive




Figure 14

100 ng/mL Mix, Positive ESI




Figure 15

ESI Negative Calibration Solution										
Name	Formula (M)	CAS	Sigma	QuantaBiodesign	FW	molarity				
7-aminoheptanoic acid	C7H15NO2	929-17-9	284637		145.2	6.89E-06				
amino-dPEG 4-acid	C11H23NO6			10244	265.3	9.42E-07				
amino-dPEG 6-acid	C15H31NO8			10067	353.41	7.07E-07				
sulfinpyrazone	C23H20N2O3S	57-96-5	S9509		404.48	1.24E-07				
amino-dPEG 8-acid	C18H39NO10	756526-04-2		10277	441.51	2.83E-07				
amino-dPEG 12-acid	C27H55NO14	756526-07-5		10287	617.72	4.61E-07				
amino-dPEG 16-acid	C35H71NO18			10297	793.93	7.87E-07				

Figure 16

1 ug/mL Mix, Negative ESI




Figure 17

ESI Negative




Figure 18

Stability of Positive Ion Solution – 5 weeks at 25°C




Figure 19

Stability of Negative Ion Solution – 5 weeks at 25°C



Figure 20

Amino-dPEG₄-acid



Mol Wt: 265.26; single compound
dPEG[®] Spacer is 15 atoms and 18.0 Å

Amino-dPEG₅-acid



Mol Wt: 353.47; single compound
dPEG[®] Spacer is 22 atoms and 25.1 Å

Amino-dPEG₆-acid



Mol Wt: 441.57; single compound
dPEG[®] Spacer is 28 atoms and 32.2 Å

Amino-dPEG₇-acid



Mol Wt: 529.72; single compound
dPEG[®] Spacer is 35 atoms and 39.5 Å

Amino-dPEG₈-acid




Mol Wt: 703.93; single compound
dPEG[®] Spacer is 51 atoms and 60.7 Å

Amino-dPEG₉-acid



Mol Wt: 878.14; single compound
dPEG[®] Spacers are 64 atoms and 75.2 Å approx

Amino-dPEG₁₀-acid



Mol Wt: 1146.35; single compound
dPEG[®] Spacer is 76 atoms and 89 Å

Amino-dPEG₁₁-acid



Mol Wt: 1476.57; single compound
dPEG[®] Spacer is 111 atoms and 132.7 Å