




(86) Date de dépôt PCT/PCT Filing Date: 2010/06/02 (87) Date publication PCT/PCT Publication Date: 2010/12/09 (45) Date de délivrance/Issue Date: 2019/07/02 (85) Entrée phase nationale/National Entry: 2011/12/05 (86) N° demande PCT/PCT Application No.: NL 2010/050332 (87) N° publication PCT/PCT Publication No.: 2010/140885 (30) Priorité/Priority: 2009/06/03 (EP09161870.2)	(51) Cl.Int./Int.Cl. G01N 33/574 (2006.01), G01N 33/58 (2006.01) (72) Inventeurs/Inventors: VAN DONGEN, JACOBUS JOHANNES MARIA, NL; ORFAO DE MATOS CORREIA E VALE, JOSE ALBERTO, NL; FLORES-MONTERO, JUAN ALEJANDRO, NL; ALMEIDA PARRA, JULIA MARIA, ES; VAN DER VELDEN, VINCENT HENRICUS JOHANNES, NL; BOETTCHER, SEBASTIAN, DE; RAWSTRON, ANDREW CRAIG, GB; ... (73) Propriétaire/Owner:
--	---

(54) **Titre : PROCÉDES, REACTIFS ET KITS D'IMMUNOPHENOTYPAGE PAR CYTOMETRIE EN FLUX**
(54) **Title: METHODS, REAGENTS AND KITS FOR FLOW CYTOMETRIC IMMUNOPHENOTYPING**



(57) **Abrégé/Abstract:**
The invention relates to the field of flow cytometry and more particularly to a panel of antibody reagents conjugated to fluorescent compounds. Provided are reagent compositions, comprising at least eight distinct fluorochrome-conjugated antibodies comprising a set of at least three identification antibodies for the identification of a leukocyte population of interest and at least four characterization antibodies for further characterization and/or classification of said leukocyte population. Also provided are kits and methods related to the reagent compositions.

- (72) **Inventeurs(suite)/Inventors(continued):** DE TUTE, RUTH MARY, GB; LHERMITTE, LUDOVIC BERNARD SIMON, FR; ASNAFI, VAHID, FR; MEJSTRIKOVA, ESTER, CZ; SZCZEPANSKI, TOMASZ, PL; MONTEIRO DA SILVA LUCIO, PAULO JORGE, PT; AYUSO, MARTA MARIN, ES; PEDREIRA, CARLOS EDUARDO, BR
- (73) **Propriétaires(suite)/Owners(continued):**ERASMUS UNIVERSITY MEDICAL CENTER ROTTERDAM, NL
- (74) **Agent:** NORTON ROSE FULBRIGHT CANADA LLP/S.E.N.C.R.L., S.R.L.

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property Organization
International Bureau(43) International Publication Date
9 December 2010 (09.12.2010)(10) International Publication Number
WO 2010/140885 A1

(51) International Patent Classification:

G01N 33/574 (2006.01) *G01N 33/58* (2006.01)

(21) International Application Number:

PCT/NL2010/050332

(22) International Filing Date:

2 June 2010 (02.06.2010)

(25) Filing Language:

English

(26) Publication Language:

English

(30) Priority Data:

09161870.2 3 June 2009 (03.06.2009) EP

(71) Applicant (for all designated States except US): **ERASMUS UNIVERSITY MEDICAL CENTER ROTTERDAM** [NL/NL]; Dr. Molewaterplein 50, NL-3015 GE Rotterdam (NL).

(72) Inventors; and

(75) Inventors/Applicants (for US only): **VAN DONGEN, Jacobus Johannes Maria** [NL/NL]; c/o Erasmus University Medical Center Rotterdam, Dr. Molewaterplein 50, NL-3015 GE Rotterdam (NL). **ORFAO DE MATOS CORREIA E VALE, José Alberto** [PT/NL]; c/o Erasmus University Medical Center Rotterdam, Dr. Molewaterplein 50, NL-3015 GE Rotterdam (NL). **FLORES-MONTERO, Juan** [ES/NL]; c/o Erasmus University Medical Center Rotterdam, Dr. Molewaterplein 50, NL-3015 GE Rotterdam (NL). **ALMEIDA PARRA, Julia Maria** [ES/ES]; C / La Paz, 4, E-37193 Cabrerizos, Salamanca (ES). **VAN DER VELDEN, Vincent Henricus Johannes** [NL/NL]; Jachthoorn 63, NL-2907 SE Capelle aan den IJssel (NL). **BÖTTCHER, Sebastian** [DE/DE]; Rosenweg 4, D-24107 Stampe (DE). **RAWSTRON, Andrew Craig** [GB/GB]; 3 The Poplars, Newton on Ouse, York Yorkshire YO30 2BL (GB). **DE TUTE, Ruth Mary** [GB/GB]; 19 Elmwood Grove, Horbury, Wakefield, West Yorkshire Yorkshire WF4 5JH(GB). **LHERMITTE, Ludovic Bernard Simon** [FR/FR]; 7, Rue George Saché, F-75014 Paris (FR). **ASNAFI, Vahid** [FR/FR]; 5, Rue Robert Branchard, F-95870 Bezons (FR). **MEJSTRÍKOVÁ, Ester** [CZ/CZ]; Nerudova 18, 11800 Prague (CZ). **SZCZEPANSKI, Tomasz** [PL/PL]; Ul. Gen. de Gaulle'a 69/5, PL-41-800 Zabrze (PL). **MONTEIRO DA SILVA LUCIO, Paulo Jorge** [PT/PT]; Rua Das Searas 31, P-2790-151 Carnaxide (PT). **AYUSO, Marta Marín** [ES/ES]; Calle Alfonso de Castro 12, 3B, E-37005 Salamanca (ES). **PEDREIRA, Carlos Eduardo** [BR/BR]; Rua Viuva Lacerda 249, BL3, 902, Rio de Janeiro (BR).(74) Agent: **HATZMANN, M.J.**; Vereenigde, Johan de Wittlaan 7, NL-2517 JR Den Haag (NL).

(81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IS, JP, KE, KG, KM, KN, KP, KR, KZ, LA, LC, LK, LR, LS, LT, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PE, PG, PH, PL, PT, RO, RS, RU, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.

(84) Designated States (unless otherwise indicated, for every kind of regional protection available): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).


Published:

— with international search report (Art. 21(3))

[Continued on next page]

(54) Title: METHODS, REAGENTS AND KITS FOR FLOW CYTOMETRIC IMMUNOPHENOTYPING

Figure 1



(57) Abstract: The invention relates to the field of flow cytometry and more particularly to a panel of antibody reagents conjugated to fluorescent compounds. Provided are reagent compositions, comprising at least eight distinct fluorochrome-conjugated antibodies comprising a set of at least three identification antibodies for the identification of a leukocyte population of interest and at least four characterization antibodies for further characterization and/or classification of said leukocyte population. Also provided are kits and methods related to the reagent compositions.

WO 2010/140885 A1

- *before the expiration of the time limit for amending the claims and to be republished in the event of receipt of amendments (Rule 48.2(h))*

Title: Methods, reagents and kits for flow cytometric immunophenotyping.

5 The invention relates to the field of flow cytometry and more particularly to a panel of antibody reagents conjugated to fluorescent compounds. They find their use in the immunophenotypic characterization of normal, reactive, regenerating and neoplastic cells in peripheral blood (PB), bone marrow (BM), pleural effusions, ascitis, cerebrospinal fluid (CSF), vitreous humor, synovial fluid, bronchoalveolar lavage,
10 urine, spleen, liver, lymph node, and other tissue samples. Flow cytometric immunophenotyping of normal, reactive, regenerating, and malignant cells (in particular leukocytes) is currently being used for many applications in medicine, including immunology, hematology and oncology. Among others, said applications include monitoring of the immune system; diagnosis and classification of primary
15 immunodeficiencies; immunophenotyping of leukemias, lymphomas and plasma cell dyscrasias; monitoring of low frequencies of malignant leukocytes as measure for treatment effectiveness; diagnosis and monitoring of clonal disorders such as paroxysmal nocturnal hemoglobinuria (PNH) and mastocytosis; evaluation of hematopoiesis or lymphopoiesis in different clinical conditions, e.g. in healthy
20 individuals, after stem cell transplantation or gene therapy (using hematopoietic precursor cells as target); and evaluation of the composition and quality of cell products to be used for therapeutic purposes.

 The conventional diagnostic process in flow cytometric immunophenotyping is typically based on the usage of panels of antibodies. At present
25 different, overlapping panels of antibodies are recommended for specific applications. Such panels are driven either by the medical indication (e.g.: screening for cytopenias, characterization of lymphocytosis), by disease (e.g.: acute leukemia diagnosis, lymphoma diagnosis) or disease status (e.g.: diagnostic classification of ALL vs. monitoring of ALL for evaluation of treatment effectiveness). Examples of
30 recommended antibody panels are the European Leukemia Net (ELN) (1), the 2006 Bethesda International Consensus (2) panels for different subtypes of hematological malignancies, the EGIL panels for lineage assignment and subclassification of acute leukemias (3), the European Myeloma Network (EMN) panels for diagnosis,

classification and monitoring of plasma cell dyscrasias (4), the ERIC (5) and Matutes (6) score panels for immunophenotypic subclassification of chronic lymphocytic leukemia and hairy cell leukemia vs. other chronic lymphoproliferative disorders.

These recommended antibody panels are mostly comparable between
5 different countries and study groups, but never fully identical. Importantly, no or limited information has been provided by the various networks and study groups on how the antibodies should preferably be combined into a panel of antibody reagents conjugated with specific fluorescence compounds. Only simple lists of antibodies have been proposed, even when diagnostic laboratories can simultaneously stain a sample
10 aliquot with 2, 3, 4, 5, 6, 7, 8 or more antibodies conjugated with different fluorescent dyes having fluorescence emissions that can be measured separately. Variations between centres around the world are due to the usage of:

- Comparable but not identical lists of antibodies and antibody clones;
- The same antibodies conjugated with different fluorochromes that show
15 different sensitivities. For example, an antibody reagent might provide a negative result if the antibody is combined with a low-sensitive fluorochrome, while a positive result could be obtained if an antibody reagent consisting of the same antibody combined with a different, more sensitive fluorochrome, is used to stain the same cells.
- 20 - Simultaneous assessment of different combinations of antibody reagents.
- Variable and frequently suboptimal or inappropriate strategies for the identification of the cell population(s) of interest, such as the CD45 gating in acute myeloblastic leukemias.

25 As a consequence, the seemingly comparable antibody panels do not result in comparable diagnosis of clinical samples. In fact, the reproducibility between diagnostic laboratories is not more than 70%. Unfortunately, diagnoses are missed because of inappropriate discrimination between normal and malignant cells or regenerating bone marrow is erroneously diagnosed as recurrence of acute leukemia.

30

In US Patent No. 5,047,321, Loken and Terstappen disclosed a procedure for the multiparameter analysis of cellular components in PB and BM. With the procedure described, these inventors could distinguish several cellular components of

PB and BM, count the number of cells within each component, and provide a differential analysis of each of them based on the combined use of the LDS-751 (Exciton) DNA-dye, the thiazol orange (TO, Molecular Probes, Inc) RNA-dye, a fluorescently labeled anti-CD45 monoclonal antibody, forward light scatter (FSC) and
5 sideward light scatter (SSC). This approach allowed the specific identification of nucleated red cells, erythrocytes, reticulocytes, platelets, lymphocytes, monocytes, neutrophilic granulocytes, basophilic granulocytes, eosinophilic granulocytes, and precursors of all nucleated hematopoietic cells. However, the described multiparameter analysis could not specifically differentiate between normal, reactive,
10 regenerating and neoplastic cells coexisting in the same sample, neither could this procedure further characterize these groups of cells. A further refinement of this patent was described later on in U.S. Patent No 6,287,791 by Terstappen and Chen, but they did not show any improvement in characterizing the different populations of leukocytes.

15 More recently, Orfao et al. described in US. Patent No. 7,332,295 a procedure for the multidimensional leukocyte differential analysis of PB, BM and other body fluids, which specifically allowed identification of dendritic cells and their subsets in addition to nucleated red cells, lymphocytes, monocytes, neutrophilic granulocytes, basophilic granulocytes, eosinophilic granulocytes and hematopoietic
20 precursors of all nucleated cells. Furthermore, in U.S. Patent No 5,538,855, Orfao et al. described a procedure that allowed a more detailed analysis of the lymphoid compartments through the simultaneous identification of up to 12 different subsets of T, B and NK-cells in PB, BM and lymph node samples. In the patent, Orfao et al. used a combined staining for the CD3, CD19, CD56 (and/or CD16), CD4 and CD8 antigens
25 in a 3-color single staining, where pairs of monoclonal antibodies conjugated with the same fluorochrome were used. Nevertheless, through this approach they could neither further characterize the identified cell subsets nor discriminate between normal and neoplastic cells; furthermore, some relevant subsets of lymphoid and non-lymphoid cells (e.g. the TCR $\gamma\delta$ + T-cells) present in a normal PB, BM or other tissues and body
30 fluids, could not be specifically detected.

None of the referred procedures allowed directly for a more detailed characterization of the cells identified, including discrimination between normal, reactive, regenerating, and neoplastic/clonal cell populations. Such characterization

and discrimination requires the use of a greater number of stainings associated with distinguishable fluorescence emissions and of flow cytometer instruments capable of measuring a higher number of different fluorescence emissions.

Many other publications show that, based on specific immunophenotypic
5 profiles, normal cells can be discriminated from their neoplastic counterparts, and that the discrimination capabilities (sensitivity for discriminating normal from clonal/malignant cells) increase with the number of fluorochrome- conjugated antibody reagents used to stain the same sample. Some reports also show that normal and reactive cells as well as maturing cells during BM regeneration, may show
10 distinct immunophenotypic profiles, but these cells could not be clearly discriminated from malignant cells. Importantly, panels of antibody combinations that have been proposed thus far can not efficiently discriminate between normal, reactive, regenerating cell populations and their neoplastic counterparts in a systematic way and, at the same time, allow for a clear-cut discrimination between different
15 diagnostic disease categories through an extended immunophenotypic characterization of the malignant cells.

With the introduction of the new generation of multi-laser flow cytometers in diagnostic laboratories, the potential space for building different panels even with the same list of antibodies, has exponentially increased due to the possibility of using
20 a larger number of different fluorescent compounds and simultaneously measure a higher number of fluorescence emissions in single cells. In U.S. Patent No. 7,321,843, Orfao, Pedreira and Sobral da Costa propose a new approach based on mathematical calculation procedures (e.g. the closest neighbour principle), where flow cytometry list mode data files containing information for each single cell measured about an
25 unlimited number of parameters could be generated, after several aliquots of a sample stained with different, partially overlapping, combinations of monoclonal antibodies had been measured in a flow cytometer. The utility of this procedure could only be maximised with the increased multicolour capabilities of the most recent flow cytometers. This is because these procedures require the design of efficient panels of
30 combinations of monoclonal antibody reagents where common backbone markers are combined with other additional markers which vary for each combination of monoclonal antibodies in a panel. In parallel, the same authors also described (U.S. Patent No. 7,507,548) a procedure for comparison of flow cytometry data from a case

against a reference data file which would allow direct comparison of e.g. normal vs. neoplastic cells. However, in both patents, the authors failed to propose a comprehensive antibody panel that could be thoroughly and systematically applied to the diagnosis, classification, staging and/or monitoring of leukaemia, lymphoma and
5 plasma cell dyscrasias.

Thus, existing protocols do not teach about how to specifically combine the widely recommended monoclonal antibodies in a panel; they propose different panels for the diagnostic classification and monitoring of minimal disease; they focus in single, uniform disease groups; they show a limited efficiency once widely tested; they
10 do not provide a way for direct combination of the information measured for several aliquots from a single sample; and/or they have failed in systematically distinguishing normal/reactive versus clonal/neoplastic cells based on their immunophenotypic properties, in specific diseases such as chronic T and NK lymphoproliferative disorders (e.g. in the differential diagnoses of large granular lymphocytic leukemia).

15 The present inventors recognized these difficulties, and set out to design improved, well-defined antibody panels that avoid misinterpretation and over-interpretation of results. After careful selection of the relevant markers, design of appropriate combinations of antibodies in multi-color tubes, and the selection of suited fluorochromes (based on need for brightness, compensation, stability, etc.), a
20 set of antibody reagents was developed. The studies were complemented with extensive multicentric evaluation of the consensus panels in order to reshape and achieve an optimal efficiency.

Accordingly, the invention provides in one embodiment a reagent composition for flow cytometric immunophenotyping of leukocytes comprising at least
25 eight distinct fluorochrome-conjugated antibodies comprising a set of at least three identification antibodies for the identification of a leukocyte population of interest and at least four characterization antibodies for further characterization and/or classification of said leukocyte population. A reagent composition as provided herein comprises a panel of antibodies directed against one of the following combinations of
30 markers: (a) CD20, CD4, CD45, CD19, Ig λ , CD8, Ig κ , CD56, TCR $\gamma\delta$, CD3 and CD38, wherein the antibody within either one of the pairs CD20/CD4, Ig λ /CD8 and CD19 /TCR $\gamma\delta$ is conjugated to the same fluorochrome (LST reagent composition); (b) CD20, CD45, CD8, Ig λ , CD56, Ig κ , CD4, CD3, CD14 and CD38, wherein the antibody within

either one of the pairs CD8/Igλ, CD56/Igκ and CD3/CD14 is conjugated to the same fluorochrome (SST reagent composition); (c) CD45, CD138, CD38, CD56, β2micro, CD19, cyIgκ and cyIgλ (PCST reagent composition) or (d) cyCD3, CD45, cyMPO, cyCD79a, CD34, CD19, CD7 and CD3 (ALOT reagent composition).

5 In a preferred embodiment, the composition comprises monoclonal antibodies. CD stands for cluster designation and is a nomenclature for the identification of specific cell surface antigens defined by monoclonal antibodies. Abbreviations used are as follows: cyIgκ = cytoplasmic IgG kappa chain; cyIgλ = cytoplasmic IgG lambda chain; β2micro = β2 microglobulin; cyMPO = cytoplasmic myeloperoxidase. (Monoclonal) antibodies against the indicated markers can be
10 commercially obtained from various companies, including Becton/Dickinson (BD) Biosciences, Dako, Beckman Coulter, CYTOGNOS, Caltag, Pharmingen, Exbio, Sanquin, Invitrogen, and the like.

Suitable fluorochromes for conjugating antibodies are known in the art. As will
15 be understood, the fluorochromes used within a reagent composition should be distinguishable by flow cytometry. The fluorochromes are preferably selected for brightness, limited spectral overlap and limited need for compensation, stability, etc. The following panel of fluorochromes is of particular use in a reagent composition according to the invention: (1) pacific blue (PacB) or Horizon V450, (2) pacific orange
20 (PacO) or AMCA, (3) fluorescein isothiocyanate (FITC) or Alexa488, (4) phycoerythrin (PE), (5) peridinin chlorophyl protein/cyanine 5.5 (PerCP-Cy5.5), PerCP or PE-TexasRed, (6) phycoerythrin/cyanine7 (PE-Cy7), (7) allophycocyanine (APC) or
Alexa647, and (8) allophycocyanine/H7 (APC-H7), APC-Cy7, Alexa680 or Alexa700.
After multiple testing rounds, the present inventors observed that very good results
25 can be obtained if the following fluorochromes are chosen: Pacific Blue or Horizon V450, Pacific Orange, fluorescein isothiocyanate (FITC), phycoerythrin (PE), peridinin chlorophyl protein/cyanine 5.5 (PerCp-Cy5.5), PE-Cy7, allophycocyanine (APC), and APC-H7.

The expression "wherein the antibody within either one of the pairs is conjugated
30 to the same fluorochrome" is meant to indicate that both antibodies of the first pair are conjugated to fluorochrome A and that both antibodies of the second pair are conjugated to fluorochrome B. Thus, within each pair the fluorochromes are the same but between different pairs the fluorochromes are distinguishable.

Each of the reagent compositions can be used as such, e.g. for the screening of a lymphoid disease. See Figure 2 for a schematic overview of exemplary applications of the various reagent compositions. The invention thus also relates to diagnostic kits comprising one or more reagent compositions. However, the compositions are also
 5 advantageously used in combination with one or more further reagent compositions, in particular reagent compositions designed for the further screening and classification of the disease. The expression "in combination with" does not refer to the physical combination or mixing of the reagent compositions, but to their application in separate (consecutive) analysis steps and combination of the data thus
 10 obtained. For example, a screening tube used in combination with a characterization tube involves two separate analytical steps on separate aliquots of the same biological sample, each using one of the reagent compositions, followed by data recording and evaluation.

Therefore, the invention also relates to a set of at least two reagent
 15 compositions, said set comprising a reagent composition as described herein above, and at least one further reagent composition comprising distinct fluorochrome-conjugated antibodies. Thus, both reagent compositions comprise a distinct panel of antibodies, although some antibodies might be present in both compositions. It is very convenient if the panel of distinct fluorochromes is essentially the same for each of the
 20 reagent compositions, and that up to eight different fluorochromes are used in total.

In one embodiment, a set of at least two reagent compositions comprises a reagent composition according to (a) and/or (b) and/or (c) defined herein above (i.e. LST reagent and/or SST reagent and/or PCST reagent), together with at least one further reagent composition comprising distinct fluorochrome-conjugated antibodies,
 25 preferably directed against one of the following combinations of markers:

- (i) CD20, CD45, CD23, CD10, CD79b, CD19, CD200 and CD43;
- (ii) CD20, CD45, CD31, LAIR1, CD11c, CD19, IgM and CD81
- (iii) CD20, CD45, CD103, CD95, CD22, CD19, CXCR5 and CD49d
- (iv) CD20, CD45, CD62L, CD39, HLADR, CD19, CD27 and CD31
- 30 (v) CD4, CD45, CD7, CD26, CD3, CD2, CD28 and CD8
- (vi) CD4, CD45, CD27, CCR7, CD3, CD45RO, CD45RA and CD8
- (vii) CD4, CD45, CD5, CD25, CD3, CD11c and CD8
- (viii) CD4, CD45, CD57, CD30, CD3, CD11c and CD8

- (ix) CD4, CD45, cyPerforin, cyGranzyme, CD3, CD16, CD94 and CD8
- (x) CD4, CD45, CD279, CD3 and CD8
- (xi) CD2, CD45, CD7, CD26, CD3, CD56, CD5 and CD19
- (xii) CD16, CD45, CD57, CD25, CD3, CD56, CD11c and CD19
- 5 (xiii) HLADR, CD45, cyPerforin, cyGranzyme, CD3, CD56, CD94 and CD19; or
- (xiv) CD45, CD138, CD38, CD28, CD27, CD19, CD117 and CD81

In one embodiment, a set comprises both the LST and SST reagents and at least one of the above further reagents. In another embodiment, a set comprises at least the
10 PCST reagent.

In another embodiment, a set of at least two reagent compositions comprising a first reagent composition according to (d) as defined herein above (ALOT reagent), and at least one further reagent composition comprising distinct fluorochrome-conjugated
15 antibodies directed against one of the following combinations of markers:

- (i) CD20, CD45, CD58, CD66c, CD34, CD19, CD10 and CD38
- (ii) IgK, CD45, cyIgμ, CD33, CD34, CD19, IgM, CD117 and Igλ, wherein the antibodies against IgM and CD117 are conjugated to the same fluorochrome
- 20 (iii) CD9, CD45, terminal deoxynucleotidyl transferase (TdT), CD13, CD34, CD19, CD22 and CD24
- (iv) CD21, CD45, CD15, CDw65, NG2, CD34, CD19, CD123 and CD81, wherein the antibodies against CD15 and CDw65 are conjugated to the same fluorochrome
- 25 (v) cyCD3, CD45, TdT, CD99, CD10, CD1a and CD3
- (vi) cyCD3, CD45, CD2, CD117, CD4, CD8, CD7 and CD3
- (vii) cyCD3, CD45, TCRγδ, TCRαβ, CD33, CD56, cyTCRβ and CD3
- (viii) cyCD3, CD45, CD44, HLADR, CD45RA, CD123 and CD3
- (ix) HLADR, CD45, CD16, CD13, CD34, CD117, CD11b and CD10
- 30 (x) HLADR, CD45, CD35, CD64, CD34, CD117, immune receptor expressed by myeloid cell 2 (IREM2) and CD14
- (xi) HLADR, CD45, CD36, CD105, CD34, CD117, CD33 and CD71

- (xii) HLADR, CD45, CD15, chondroitin sulfate proteoglycan (NG2), CD34, CD117, CD7 and CD38
- (xiii) HLADR, CD45, CD42a, CD61, CD203c, CD34, CD117, CD123 and CD4, wherein the antibodies against CD42a and CD61 are conjugated to the same fluorochrome
- (xiv) HLADR, CD45, CD41, CD25, CD34, CD117, CD42b and CD9
- (xv) HLADR, CD45, CD41, CD25, CD34, CD117, CD42b and CD9; or
- (xvi) HLADR, CD45, TdT, CD56, CD34, CD117, CD22 and CD19.

Exemplary sets of reagent compositions include at least one reagent composition as recited in Tables 1 and/or 6 herein below, together with at least one reagent composition as recited in any one of Tables 2 -5 and 7-9.

A further aspect relates to a diagnostic kit for flow cytometric immunophenotyping of leukocytes, wherein the kit comprises one or more sets of at least two reagent compositions described above. The kit may in addition comprise other useful components, such as instructions for use, sample preparation reagent, buffer, and/or control samples.

Reagent compositions, sets and diagnostic kits provided herein find their application in various fields. For example, the proposed panels can be applied as a whole or only partially depending on the nature of the sample, medical indication or the specific goal. The panel may use specific monoclonal antibody reagents by a single or several different manufacturers, in combination with different cell preparation techniques for staining of cell surface only and/or cell surface plus intracellular markers. Similarly, the reagent compositions may be upgraded into a panel including combinations containing monoclonal antibody reagents conjugated with > 8 different fluorochromes where the same backbone markers are maintained and combined with additional or similar markers. An antibody composition may also act as a core panel to which new combinations are added with different goals, including monitoring of the immune system.

In one embodiment, there is provided a diagnostic kit for the identification and characterization of mature lymphoid cells, comprising the SST and/or LST and/or PCST reagent composition together with at least one reagent composition for

detecting B-cell chronic lymphoproliferative disorder (B-CLPD) comprising antibodies against CD20, CD45, CD23, CD10, CD79b, CD19, CD200 and CD43; CD20, CD45, CD31, LAIR1, CD11c, CD19, IgM and CD81; CD20, CD45, CD103, CD95, CD22, CD19, CXCR5 and CD49d; or CD20, CD45, CD62L, CD39, HLADR, CD19, CD27 and
 5 CD31.

For example, there is provided a diagnostic kit for the identification and characterization of mature lymphoid cells, comprising the SST and/or LST and/or the PCST reagent composition together with at least one reagent composition for detecting T-cell chronic lymphoproliferative disorder (T-CLPD) comprising antibodies
 10 against CD4, CD45, CD7, CD26, CD3, CD2, CD28 and CD8; CD4, CD45, CD27, CCR7, CD3, CD45RO, CD45RA and CD8; CD4, CD45, CD5, CD25, CD3, CD11c and CD8; CD4, CD45, CD57, CD30, CD3, CD11c and CD8; CD4, CD45, cyPerforin, cyGranzyme, CD3, CD16, CD94 and CD8; or CD4, CD45, CD279, CD3 and CD8.

Another exemplary diagnostic kit for the identification and characterization of mature lymphoid cells comprises the SST and/or LST and/or PCST reagent composition together with at least one reagent composition for detecting NK-cell chronic lymphoproliferative disorder (NK-CLPD) comprising antibodies against CD2, CD45, CD7, CD26, CD3, CD56, CD5 and CD19; CD16, CD45, CD57, CD25, CD3, CD56,
 20 CD11c and CD19; HLADR, CD45, cyPerforin, cyGranzyme, CD3, CD56, CD94 and CD19; CD45, CD138, CD38, CD56, β 2micro, CD19, cyIgk and cyIg λ ; or CD45, CD138, CD38, CD28, CD27, CD19, CD117 and CD81.

In yet another embodiment, there is provided a diagnostic kit for the identification and characterization of mature lymphoid cells, comprising the SST
 25 and/or LST and/or PCST reagent composition together with at least one reagent composition for detecting plasma cell dyscrasias (PCD) comprising antibodies against CD45, CD138, CD38, CD56, β 2micro, CD19, cyIgk and cyIg λ ; or CD45, CD138, CD38, CD28, CD27, CD19, CD117 and CD81.

Also provided is a diagnostic kit for identification and characterization of
 30 immature lymphoid cells, comprising the ALOT reagent composition together with at least one reagent composition for detecting B-cell precursor ALL (BCP-ALL) comprising antibodies against CD20, CD45, CD58, CD66c, CD34, CD19, CD10 and CD38; Igk, CD45, cyIg μ , CD33, CD34, CD19, IgM, CD117 and Ig λ , wherein the

antibodies against IgM and CD117 are conjugated to the same fluorochrome; CD9, CD45, TdT, CD13, CD34, CD19, CD22 and CD24; or CD21, CD45, CD15, CDw65, NG2, CD34, CD19, CD123 and CD81, wherein the antibodies against CD15 and CDw65 are conjugated to the same fluorochrome.

5 In a further aspect, a diagnostic kit for identification and characterization of immature lymphoid cells comprises the ALOT reagent composition together with at least one reagent composition for detecting T-cell precursor ALL (T-ALL) comprising antibodies against cyCD3, CD45, TdT, CD99, CD10, CD1a and CD3; cyCD3, CD45, CD2, CD117, CD4, CD8, CD7 and CD3; cyCD3, CD45, TCR $\gamma\delta$, TCR $\alpha\beta$, CD33, CD56,
10 TCR β F1 and CD3; or cyCD3, CD45, CD44, HLADR, CD45RA, CD123 and CD3.

In still a further aspect, a diagnostic kit is provided for identification and characterization of myeloid cells, comprising the ALOT reagent composition together with at least one reagent composition for detecting acute myeloid leukemia (AML), myelodysplastic syndrome (MDS)/ chronic myeloproliferative disorder (MPD),
15 comprising antibodies against HLADR, CD45, CD16, CD13, CD34, CD117, CD11b and CD10; HLADR, CD45, CD35, CD64, CD34, CD117, IREM2 and CD14; HLADR, CD45, CD36, CD105, CD34, CD117, CD33 and CD71; HLADR, CD45, CD15, NG2, CD34, CD117, CD7 and CD38; HLADR, CD45, CD42a, CD61, CD203c, CD34, CD117, CD123 and CD4, wherein the antibodies against CD42a and CD61 are conjugated to
20 the same fluorochrome; HLADR, CD45, CD41, CD25, CD34, CD117, CD42b and CD9; HLADR, CD45, CD41, CD25, CD34, CD117, CD42b and CD9; or HLADR, CD45, TdT, CD56, CD34, CD117, CD22 and CD19.

The invention also relates to a method for flow cytometric immunophenotyping of leukocytes, comprising the steps of providing a biological sample comprising
25 leukocytes and contacting at least a portion (aliquot) of the sample with a reagent composition provided herein. Any type of (human) sample known or suspected to contain leukocytes may be used. For example, the sample is peripheral blood, bone marrow, tissue sample such as lymph nodes, adenoid, spleen, or liver, or other type of body fluid such as cerebrospinal fluid, vitreous fluid, synovial fluid, pleural effusions
30 or ascitis. Preferably, the method comprises contacting a first aliquot of said sample with a first reagent composition of a set according to the invention and contacting at least a second aliquot of said sample with a further reagent composition of said set; analyzing leukocytes in said aliquots in a flow cytometer; and storing and evaluating

the data obtained. Typically, step (c) comprises the use of software for data integration and multidimensional analysis of flow cytometry files. Very suitable for use in a method of the invention is the software commercially available from CYTOGNOS SL (Salamanca, Spain) under the tradename INFINICYT™. The

5 INFINICYT™ software can automatically combine the immunophenotypic information of the selected cell populations from multiple tubes according to the so-called nearest neighbor calculations in which individual cells from one aliquot of a sample are matched with corresponding individual cells from another aliquot of the same sample, according to their backbone markers and scatter profile. The

10 INFINICYT procedure can transform the herein presented 8-color EuroFlow panels into 12, 16, or ≥ 20 -color immunostainings, dependent on the number of tubes per panel and the number of backbone markers per tube. The antibody panels and the INFINICYT software can be used in combination with all currently available flow cytometers that allow 8-color immunostainings, such as FACSCanto™ II, FACS Aria,

15 LSRII, DAKO CyAn™, Gallio, etc.

LEGEND TO THE FIGURES

Figure 1. Composition of three Categories of antibody panels

20 Figure 2. Diagnostic flow diagram showing potential applications of the EuroFlow antibody panels.

DETAILED DESCRIPTION

25 1. Introduction

The present study was performed by the The EuroFlow Consortium (LSHB-CT-2006-018708) who initiated the project “Flow Cytometry for Fast and Sensitive Diagnosis and Follow-up of Haematological Malignancies”, which includes the design of standardized multicolor immunophenotyping protocols for diagnosis, classification,

30 and monitoring of leukemias, lymphomas and plasma cell dyscrasias. A key innovative component of these protocols are the EuroFlow panels of antibody combinations (EuroFlow panels). This EuroFlow invention relates among others to panels of combinations of antibody reagents (reagent compositions), which can be

used to define normal, reactive, regenerating and malignant hematopoietic cells in a standardized way. As such, they allow for a comprehensive immunophenotypic diagnosis, classification, staging and monitoring of both chronic and acute leukemias, myelodysplastic syndromes, myeloproliferative disorders, mastocytosis, paroxysmal nocturnal haemoglobinuria, lymphomas and plasma cell dyscrasias. For the first time, the proposed panels are not based on subjective expert opinions, but they have been tested prospectively and modified for improving the answer to the most frequent medical indications of flow cytometry immunophenotyping. In turn, these panels are designed in an innovative way to be applied in combination with both conventional data analysis approaches and new interactive and semi-automated data analysis procedures in which information on single cells is combined for all parameters derived from the measurement of staining a sample with the antibody panel.

In order to establish these panels, the following sequential steps were performed:

1. Evaluation of all relevant markers used in the field for their usefulness or added value
 - Markers proposed by European networks, such as ELN, EMN, ERIC;
 - Markers proposed by the 2006 Bethesda International Consensus in US;
 - New markers for mature B-cell malignancies as proposed by Rawstron et al.
2. Design and selection of combinations of markers (≥ 6), which can recognize normal vs. reactive vs. regenerating vs. abnormal/malignant cells within a specific compartment of hematopoietic cells (cell lineage, differentiation pathway, maturation stage and/or functional subset).
3. Evaluation of the proposed panels of combinations of antibody reagents in primary samples from healthy subjects and patients.
4. Repeated testing and optimization of the combinations of antibody reagents (choice of marker, choice of antibody clone, and choice of fluorochrome) based on uncommitted goals and room for improvement.
5. Evaluation of the optimized versions of the EuroFlow panels on large series of well-defined patient samples and samples from healthy controls.

Compared with pre-existing knowledge, this is the first design of a comprehensive panel (after prospectively evaluated in a multicentric way), which allows both

discrimination between normal, reactive, regenerating and clonal/neoplastic cells and classification, staging and monitoring of clonal/neoplastic haematopoietic disorders, providing a clear indication about: 1) those markers required to be stained in common for appropriate and reproducible identification of all cell populations of interest in all
5 stained aliquots of a sample, and 2) how they should be combined with further characterization markers in specific combinations of fluorochrome-conjugated antibody reagents. In addition, information about the goals of each combination is also given as indication about when and how to apply it. The invention was made only after extensive antibody panel testing and several redesigning cycles.

10

The EuroFlow panels of reagents are composed of subsets of one or multiple combinations (named tubes) of antibodies conjugated with eight or more fluorescent compounds, each of said combinations of reagents having different goals. The EuroFlow panels consist of three different categories: Category 1, Category 2 and
15 Category 3 antibody panels (Figure 1).

2. Categories of antibody panels

The following purposes are pursued by the EuroFlow panels of antibody combinations:
20

The Category 1 antibody panels aim at identification and characterization of different subsets of mature lymphoid cells, including normal, reactive, regenerating and neoplastic B-, T-, NK-cells and plasma cells, particularly in samples where a clonal and/or neoplastic lymphoid disorder is suspected because of e.g. lymphocytosis, lymph
25 node enlargement, splenomegaly, monoclonal serum component, unexplained neurological symptoms, etc. The Category 2 antibody panels aim at identification and characterization of normal, reactive, regenerating and neoplastic immature (or early maturing) T- and/or B-lymphoid cells, particularly in samples suspected of containing neoplastic lymphoid precursors. The Category 3 antibody panels aim at identification
30 and characterization of normal, reactive, regenerating and neoplastic immature, maturing and matured myeloid cells, particularly in samples suspected of containing neoplastic myeloid cells or neoplastic cells expressing myeloid-associated markers.

Within each of said three Categories, some of the combinations of antibody reagents are aimed to be used in a single-tube screening step with more broad aims, while others are more likely applied in multi-tube classification steps when more specific target populations have already been identified in the screening step.

5

2.1 Category 1 antibody panels

The Category 1 antibody panels are composed of three screening tubes aimed at the initial identification and characterization of the specific subgroups of mature lymphoid cells present in samples which contain normal or high cell counts (e.g. normal peripheral blood) and low cell counts (e.g. vitreous humor), respectively, and four different sets of multi-tube antibody combinations aiming at further characterization of B-, T-, NK-cells and plasma cells. The Category 1 screening tubes may be used to screen for the presence of clonal and neoplastic T, B- or NK-cells and plasma cells in samples with relatively high and low cell counts, respectively. Typical examples of low cell count samples are fine needle aspirates (FNA), cerebrospinal fluid (CSF), and vitreous humor.

The above mentioned screening tubes are now termed as lymphoid screening tube for high cell count samples (abbreviated as LST), small sample tube (SST) for low cell count samples, and a separate plasma cell screening tube (PCST). In turn, the other four sets of tubes are devoted, among other uses, to further characterize the identified clonal or neoplastic B-, T-, NK-cells and plasma cells in patients with different B-, T- and NK-cell chronic lymphoproliferative disorders (abbreviated as BCLPD, TCLPD and NKCLPD, respectively) and plasma cell dyscrasias (abbreviated as PCD), respectively.

25

2.2 Category 2 antibody panels

The Category 2 antibody panels consist of a screening tube devoted to the initial identification and characterization of immature lymphoid vs. non-lymphoid precursors, and two sets of tubes for further detailed characterization of B-cell and T-cell precursors. Among other applications, the screening tube may be used for the classification of acute leukemias into lymphoid vs. non-lymphoid versus undifferentiated versus biphenotypic or bilineage acute leukemias. This tube might also be used for the discrimination between normal/regenerating and

30

neoplastic/malignant B-cell precursors in the bone marrow, peripheral blood and other tissues and T-cell precursors in the thymus, respectively. Consequently, this screening tube is now termed as acute leukemia orientation tube (abbreviated ALOT).

The other two sets of tubes may be used, among other applications, for the

- 5 characterization of lymphoid blast cells showing B-cell and T-cell commitment in patients with B-cell precursor acute lymphoblastic leukemia (BCP-ALL) and T-ALL as well as in patients with acute leukemia showing an undifferentiated or biphenotypic/bilineage phenotype.

10 *2.3 Category 3 antibody panels*

Finally, the Category 3 antibody panels are composed of a screening tube which aims at discriminating between lymphoid and non-lymphoid precursor cells and which is identical to the screening tube described above in Category 2 (ALOT tube), and one additional set of tubes aiming at the characterization of maturing non-lymphoid cells

- 15 from the earliest hematopoietic precursors to mature myeloid cells, including the different maturation stages of cells committed to the erythroid, megakaryocytic, monocytic, neutrophilic, eosinophilic, basophilic, mast cell and dendritic cell lineages as well as the more immature, uncommitted hematopoietic precursors and stromal cells in case they are present in the sample (e.g. in BM). Among other uses, this set of
- 20 tubes may be applied for the characterization of neoplastic/clonal diseases and other disorders in which the cells of interest (e.g.: clonal/neoplastic, altered/aberrant cells) display myeloid differentiation such as in patients with acute myeloid leukemia (AML), myelodysplastic syndrome (MDS), chronic myeloproliferative disorder (MPD), paroxysmal nocturnal haemoglobinuria (PNH), mastocytosis, idiopathic cytopenia of
- 25 undetermined significance (ICUS)- and/or acute leukemia with myeloid-associated markers as it occurs in biphenotypic and bilineage acute leukemias. This group of tubes is abbreviated from now on as AML/MDS/MPD.

2.4 Types of markers used in EuroFlow panels

- 30 In each tube of the EuroFlow panels, two types of antibody reagents are combined: 1) reagents which mainly aim at identifying the precise cell populations of interest present in the sample (Backbone markers), which also provide additional information about the phenotypic characteristics of said cell populations, and; 2) reagents mainly

devoted to further characterization/classification of said cell populations as well as other groups of cells in the sample (Characterization markers).

Typically, backbone markers are repeated in each tube of the following sets of tubes: AML/MDS/MPD, BCP-ALL, T-ALL, B-CLPD, T-CLPD, NK-CLPD, PCD). In
 5 addition, backbone markers of the B-CLPD set of tubes are also common to the LST and SST screening tubes, backbone markers of the PCD set of tubes are common to the PCST, and backbone markers of the BCP-ALL and the T-ALL sets of tubes are also common to the ALOT; finally, two of the backbone markers in the AML panel (i.e.: CD34 and CD45) are also common to the ALOT tube.

10 In turn, characterization markers are combined in a comprehensive way in each tube, so that they would allow distinguishing normal, reactive regenerating vs clonal/neoplastic cells, even when present in low numbers and in case of neoplastic cells they would allow further diagnosis, subclassification, staging and monitoring of acute and chronic leukemias, lymphomas and plasma cell dyscrasias. For this purpose
 15 in multi-tube panels, each tube is devoted to a specific aim related to full characterization and monitoring of a disease entity if combined with the information from the corresponding screening tube (e.g. for the diagnosis, staging and monitoring of CLL, tube#4 of the B-CLPD multi-tube panel in combination with the LST will suffice) or for specific disease –associated information (e.g. tube# 24 in the BCP-ALL
 20 multi-tube panel shown in Table 7 is specifically devoted for distinguishing between normal/reactive and regenerating versus malignant B-cell precursors and minimal residual disease monitoring in BCP-ALL).

2.5 INFINICYT software tools for optimal application of the EuroFlow panels

25 The antibody panels according to the invention can be used in combination with the INFINICYT software tools which are commercially available. The software is based on recently described procedures for generating files with an unlimited number of parameters through merging data files and calculating the information derived from the measurement of markers in one sample aliquot to the individual cells
 30 measured in other aliquots of the same sample, using different combinations of antibody reagents which only have partial overlap (US 7,321,843), as well as for comparisons between different samples or different groups of samples (US 7,507,548).

3. Composition of the Category 1 antibody panels for mature lymphoid cells

5 In the multi-tube combinations of antibody reagents of the Category 1 antibody panels for mature lymphoid cells the selected common backbone markers vary between the four different panels and consist of (antibody CD plus fluorochrome compound number): 1) CD20-1, CD45-2 and CD19-6 for B-CLPD; 2) CD4-1, CD45-2, CD3-5 and CD8-8 for T-CLPD; 3) CD45-2, CD3-5, CD56-6 and CD19-8 for NK-CLPD, and; 4) 10 CD45-1, CD138-2, CD38-3 and CD19-6 for PCD. In addition, the backbone markers used in the B-CLPD tube are also used in the LST and SST screening tubes, and the backbone markers of the PCD tubes are also used in the PCST.

The backbone markers in each of these multi-tube combinations of antibody reagents 15 aim at providing delineation of the groups of cells of interest and the specific identification of neoplastic cells in diagnostic and follow-up samples containing enough numbers of tumor cells. In addition, for each of the four multi-tube combinations, backbone markers are combined with a variable number of additional characterization markers which further contribute to the discrimination between 20 normal, reactive, regenerating and neoplastic/clonal B-, T, NK-cells and plasma cells, respectively, even when these are present in minimal numbers (e.g. minimal residual disease monitoring and disease staging), as well as for the distinction between clonal/neoplastic cells from different disease categories.

25 *3.1. The lymphoid screening tube (LST), the small sample tube (SST), and the plasma cell screening tube (PCST)*

The EuroFlow LST reagent composition was designed and approved for evaluation of several suspected clinical conditions, such as lymphocytosis, lymph node enlargement, 30 splenomegaly, monoclonal serum components, unexplained neurological symptoms, unexplained cytopenias, etc. (Figure 2). The composition of an exemplary LST tube is provided in Table 1. This tube detects aberrant mature lymphocyte populations of B, T and NK lineage. However, this 8-color tube does not allow the precise diagnosis and

classification of the detected aberrant lymphocyte populations. This typically needs further characterization with the B-CLPD, T-CLPD, NK-CLPD and/or PCD tubes (see Sections 3.2, 3.3 and 3.4).

- 5 The EuroFlow SST reagent composition is a modified version of EuroFlow LST, specially designed for evaluation of small samples and samples with (very) low cell counts, such as fine needle aspirates (FNA), cerebrospinal fluid (CSF), vitreous humor, etc. For this special aim, the tube allows the unequivocal recognition of normal leukocytes present in these samples, e.g. B, T, NK cells and monocytes as well
10 as any coexisting aberrant cell population. The composition of an exemplary SST reagent is provided in Table 1.

The EuroFlow PCST reagent composition is specially designed for screening of plasma cells in order to detect aberrancies or clonality. In case of aberrant or clonal plasma
15 cells, complementary phenotypic characterization is achieved via the 2-tube PCD antibody panel.

In any of the Tables of the present application, fluorochrome number 1 corresponds to pacific blue (PacB) or Horizon V450, number 2 to pacific orange (PacO) or AMCA,
20 number 3 to fluorescein isothiocyanate (FITC) or Alexa488, number 4 to phycoerythrin (PE), number 5 to peridinin chlorophyll protein/cyanine 5.5 (PerCP-Cy5.5), PerCP or PE-Texas Red, number 6 to phycoerythrin/cyanine7 (PE-Cy7), number 7 to allophycocyanine (APC) or Alexa647, and number 8 to allophycocyanine/H7 (APC-H7), APC-Cy, Alexa680 or Alexa700. A preferred
25 combination is pacific blue (PacB), pacific orange (PacO), fluorescein isothiocyanate (FITC), phycoerythrin (PE), peridinin chlorophyll protein/cyanine 5.5 (PerCP-Cy5.5), phycoerythrin/cyanine7 (PE-Cy7), allophycocyanine (APC), and allophycocyanine/H7 (APC-H7).

30

35

TABLE 1. Single tube EuroFlow screening combinations for normal and malignant/ aberrant mature lymphoid cells

5

Fluorochrome	1	2	3	4	5	6	7	8
Tube#								
1 (LST)	CD20 ^{BB} and CD4	CD45 ^{BB}	Igλ and CD8	Igκ and CD56	CD5	CD19 ^{BB} and TCRγδ	CD3	CD38
2 (SST)	CD20 ^{BB}	CD45 ^{BB}	CD8 and Igλ	CD56 and Igκ	CD4	CD19 ^{BB}	CD3 and CD14	CD38
3 (PCST)	CD45 ^{BB}	CD138 ^{BB}	CD38 ^{BB}	CD56	β2micro	CD19 ^{BB}	CyIgκ	CyIgλ

BB =These markers act as Backbone markers when LST or SST results are combined with the results obtained with multi-tube EuroFlow classification combination for B-CLPD (see Table 2) or when the PCST results are combined with the results obtained with the 2-tube EuroFlow classification combination for PCD (see Table 5); LST = lymphoid screening tube; SST = small sample tube; PCST= plasma cell screening tube.

3.2. The multi-tube antibody panel for B-cell chronic lymphoproliferative disorders (B-CLPD)

The B-CLPD is designed to classify mature B-cell malignancies according to WHO entities based on flow cytometric data only (see Table 2). Information obtained using the LST simultaneously or sequentially has to be integrated into the B-CLPB panel (e.g. via INFINICYT software). The B-CLPD panel is designed to work in cases in which the malignant B cell population can be purified to > 90% using the backbone markers CD20, CD19, and CD45, regardless of the cell material analyzed.

The panel is designed to work in a modular way, i.e. it is not necessary to stain the whole panel if the pre-test probability for a particular B-cell malignancy is high. In those instances the panel will allow to diagnose a particular entity using a reduced number of tubes. For example, the LST tube#1 plus tube# 5 are sufficient to diagnose

5 CLL with a very high positive predictive value (PPV).

10 **TABLE 2. Multi-tube EuroFlow classification combinations for B-cell chronic lymphoproliferative disorders (B-CLPD)**

	1	2	3	4	5	6	7	8	Aim
Tube#									
4	CD20 ^{BB} and CD4	CD45 ^{BB}	Igλ and CD8	Igκ and CD56	CD5	CD19 ^{BB} and TCRγδ	CD3	CD38	LST: Detection of (almost) all mature B-cell malignancies
5	CD20 ^{BB}	CD45 ^{BB}	CD23	CD10	CD79b	CD19 ^{BB}	CD200	CD43	Identification, of all CLL cases
6	CD20 ^{BB}	CD45 ^{BB}	CD31	LAIR1	CD11c	CD19 ^{BB}	IgM	CD81	Identification of all HCL cases; characterization of benign B-, T-, NK-cells
7	CD20 ^{BB}	CD45 ^{BB}	CD103	CD95	CD22	CD19 ^{BB}	CXCR5	CD49d	Identification of DLBCL, FL, MZL, LPL
8	CD20 ^{BB}	CD45 ^{BB}	CD62L	CD39	HLADR	CD19 ^{BB}	CD27	CD31	Identification of DLBCL, FL, MZL, LPL

BB = Backbone marker; Tube 4 is identical to the LST (see Table 1). The described tubes can also be successfully applied for disease staging and monitoring.

15

3.3. The multi-tube antibody panel for T-cell chronic lymphoproliferative disorders (T-CLPD)

20 The EuroFlow T-CLPD aims for diagnosis and classification of mature T-cell malignancies. Also for this panel, information obtained with the LST simultaneously and sequentially is preferably integrated with T-CLPD tubes (e.g, via harmonization with the INFINICYT software). The panel is designed to work in cases in which the malignant T cell population can be purified to > 90% using the backbone-markers

25 CD3, CD4, CD8, and CD45, regardless of the cell material analyzed. The combination

of LST and T-CLPD tubes can detect T-cell malignancies of both TCR $\alpha\beta$ and TCR $\gamma\delta$ lineages.

5

TABLE 3. Multi-tube EuroFlow classification combinations for T- cell chronic lymphoproliferative disorders (T-CLPD)

Tube#	1	2	3	4	5	6	7	8	Aim
9	CD4 ^{BB}	CD45 ^{BB}	CD7	CD26	CD3 ^{BB}	CD2	CD28	CD8 ^{BB}	Phenotypic characterization; Identification of Sezary syndrome
10	CD4 ^{BB}	CD45 ^{BB}	CD27	CCR7	CD3 ^{BB}	CD45RO	CD45RA	CD8 ^{BB}	Phenotypic characterization; Assessment of maturation stage
11	CD4 ^{BB}	CD45 ^{BB}	CD5	CD25	CD3 ^{BB}	HLADR	cyTCL1	CD8 ^{BB}	Phenotypic characterization; Identification of T-PLL
12	CD4 ^{BB}	CD45 ^{BB}	CD57	CD30	CD3 ^{BB}		CD11c	CD8 ^{BB}	Phenotypic characterization; Cytotoxic phenotype and identification of anaplastic lymphoma
13	CD4 ^{BB}	CD45 ^{BB}	cyPerforin	cyGranzyme	CD3 ^{BB}	CD16	CD94	CD8 ^{BB}	Phenotypic characterization; Assessment of cytotoxic-associated phenotypes; Identification of T-LGL
14	CD4 ^{BB}	CD45 ^{BB}		CD279	CD3 ^{BB}			CD8 ^{BB}	Identification of lymphomas derived from follicular helper T cells (angioimmunoblastic T cell lymphomas)

10

BB = Backbone marker; the tubes can also be successfully applied for disease staging and monitoring.

15

3.4. The multi-tube antibody panel for NK-cell chronic lymphoproliferative disorders (NK-CLPD)

The EuroFlow NK-CLPD tube aims at the discrimination between aberrant and normal/reactive NK-cells. The NK-CLPD panel uses four Backbone markers: CD45-2, CD3-5, CD56-6 and CD19-8. (NK-CLPD)

5 TABLE 4. Multi-tube EuroFlow classification combinations for NK-cell chronic lymphoproliferative disorders (NK-CLPD)

Tube#	1	2	3	4	5	6	7	8	Aim; MRD
15	CD2	CD45 ^{BB}	CD7	CD26	CD3 ^{BB}	CD56 ^{BB}	CD5	CD19 ^{BB}	Detection of aberrant NK cell phenotype
16	CD16	CD45 ^{BB}	CD57	CD25	CD3 ^{BB}	CD56 ^{BB}	CD11c	CD19 ^{BB}	Detection of aberrant NK cell phenotype
17	HLADR	CD45 ^{BB}	cyPerforin	cyGranzyme	CD3 ^{BB}	CD56 ^{BB}	CD94	CD19 ^{BB}	Detection of aberrant NK cell phenotype; Assessment of cytotoxic effector phenotype

10 BB = Backbone marker; the tubes can also be successfully applied for disease staging and monitoring.

15 *3.5. The multi-tube antibody panel for plasma cell dyscrasias (PCD)*

The EuroFlow PCD panel comprises two tubes with four Backbone markers: CD45-1, CD138-2, CD38-3 and CD19-6 for PCD. (Table 5); tube #18 is identical to the PCST (tube #3) in Table 1. The PCD panel aims at the identification and enumeration of plasma cells as well as at the discrimination between normal polyclonal plasma cells such as in reactive plasmacytosis versus aberrant monoclonal plasma cells such as in monoclonal gammopathies of undetermined significance (MGUS), smoldering and symptomatic multiple myeloma (MM), plasma cell leukemias (PCL), amyloidosis, and extramedullary plasmacytoma, and the differential diagnosis among these plasma cell dyscrasias. In combination with the EuroFlow LST and B-CLPD panels, this multi-tube antibody panel will also contribute to the diagnosis of other plasma cell dyscrasias such as Waldenström's macroglobulinemia and lymphoplasmacytic lymphoma (LPL).

5 TABLE 5. Multi-tube EuroFlow classification combinations for plasma cell dyscrasias (PCD)

Tube#	1	2	3	4	5	6	7	8	Aim;
18	CD45 ^{BB}	CD138 ^{BB}	CD38 ^{BB}	CD56	β2micro	CD19 ^{BB}	cyIgk	cyIgλ	PCST: Detection of aberrant and clonal plasma cells
19	CD45 ^{BB}	CD138 ^{BB}	CD38 ^{BB}	CD28	CD27	CD19 ^{BB}	CD117	CD81	Complementary phenotypic characterization and evaluation of markers with potential prognostic impact

BB = Backbone marker. The described tubes can also be successfully applied for disease staging and monitoring; Tube #18 is identical to the PCST (Tube #3 in Table 1).

4. Composition of the Category 2 panel for immature lymphoid cells

15

4.1. The acute leukemia orientation tube (ALOT)

The EuroFlow ALOT tube was designed for assessment of the nature of immature blast cell populations in acute leukemia patients (B, T versus non-lymphoid, undifferentiated or mixed phenotype acute leukemias) and consequent orientation towards the most appropriate complementary antibody panel(s): BCP-ALL, T-ALL, and/or AML/MDS/MPD. The composition of the ALOT tube is provided in Table 6; the markers CD45-2, CD34-5 and CD19-6 act as backbone markers when the information of the ALOT is combined with the BCP-ALL panel (e.g.: using the INFINICYT software) and the cyCD3-1, CD45-2 and CD3-8 markers act also as common backbone markers when the information of the ALOT is combined with the T-ALL panel; in addition CD45-2 and CD34-5 are also common backbone markers to the AML/MDS/MPD multi-tube panel.

5 **TABLE 6. Single tube EuroFlow combination for lineage assignment of blast cells in acute leukemias**

Tube#	1	2	3	4	5	6	7	8
20 (ALOT)	cyCD3 ^{BB}	CD45 ^{BB}	cyMPO	cyCD79a	CD34 ^{BB}	CD19 ^{BB}	CD7	CD3 ^{BB}

10 BB =These markers act as Backbone markers when results are combined with the results obtained with multi-tube EuroFlow classification combinations for BCP-ALL or T-ALL (see Table 7 and Table 8, respectively);

15 The ALOT is not suitable for exclusion of a hematological malignancy, because the ALOT antibody combination is not sufficient for that purpose. However, when the ALOT is combined with the LST and 4 tubes of the AML/MDS protocol (tube # 29, 30, 31 and 32), virtually all types of hematological malignancies can be detected (not classified) or excluded (Figure 2).

4.2. *The multi-tube antibody panel for B-cell precursor acute lymphoblastic leukemia (BCP-ALL)*

20 The EuroFlow BCP-ALL tube aims at the recognition and classification of all classically defined BCP-ALL (pro-B-ALL, common-ALL, pre-B-ALL) or alternative BCP-ALL classifications, including immunophenotypic classifications associated with well-defined molecular aberrations, such as specific fusion genes. The information
25 obtained with the BCP-ALL tube set needs to be combined with ALOT, based on the backbone markers CD45-2, CD34-5 and CD19-6 and using appropriate data analysis software (c.g. via the INFINICYT software) (Table 7).

30

5 **TABLE 7. Multi-tube EuroFlow classification combinations for B-cell precursor ALL (BCP-ALL)**

Tube#	1	2	3	4	5	6	7	8	Aim
21	CD20	CD45 ^{BB}	CD58	CD66c	CD34 ^{BB}	CD19 ^{BB}	CD10	CD38	Diagnosis and classification of BCP-ALL; Detection of LAP markers
22	Igκ	CD45 ^{BB}	CyIgμ	CD33	CD34 ^{BB}	CD19 ^{BB}	IgM and CD117	Igλ	Diagnosis and classification of BCP-ALL; Detection of phenotypes associated with molecular aberrations
23	CD9	CD45 ^{BB}	TdT	CD13	CD34 ^{BB}	CD19 ^{BB}	CD22	CD24	Diagnosis and classification of BCP-ALL; Detection of phenotypes associated with molecular aberrations; Detection of LAP markers
24	CD21	CD45 ^{BB}	CD15 and CDw65	NG2	CD34 ^{BB}	CD19 ^{BB}	CD123	CD81	Subclassification of BCP-ALL; Detection of LAP markers

BB = Backbone markers. The tubes can also be successfully applied for disease staging and monitoring.

10

4.3. *The multi-tube antibody panel for T-cell acute lymphoblastic leukemia (T-ALL)*

15 The EuroFlow T-ALL panel consists of four tubes and uses cyCD3-1, CD45-2, and CD3-8 as backbone markers (see Table 8) in common with the ALOT. The T-ALL panel aims at the recognition and classification of all classically defined T-ALL (immature T-ALL, common thymocytic T-ALL, mature T-ALL) or alternative T-ALL classification, e.g. based on TCR protein expression (cyTCRβ, TCRαβ, TCRγδ) or based
20 on association with well-defined molecular aberrations.

5

TABLE 8. Multi-tube EuroFlow classification combinations for T-ALL

Tube#	1	2	3	4	5	6	7	8	Aim
25	cyCD3 ^{BB}	CD45 ^{BB}	TdT	CD99	CD5	CD10	CD1a	CD3 ^{BB}	Diagnosis of T-ALL; Identification of LAP markers
26	cyCD3 ^{BB}	CD45 ^{BB}	CD2	CD117	CD4	CD8	CD7	CD3 ^{BB}	Diagnosis of T-ALL, classification of T-ALL, and determine the maturation stage of arrest; Identification of LAP markers
27	cyCD3 ^{BB}	CD45 ^{BB}	TCR $\gamma\delta$	TCR $\alpha\beta$	CD33	CD56	Cy TCR β	CD3 ^{BB}	Diagnosis of T-ALL and determine the maturation stage of arrest; Identification of LAP markers
28	cyCD3 ^{BB}	CD45 ^{BB}	CD44	CD13	HLADR	CD45RA	CD123	CD3 ^{BB}	Subclassification of T-ALL

10 BB = Backbone markers. The tubes can also be successfully applied for disease staging and monitoring.

5. Composition of the Category 3 panel for myeloid cells

15

The EuroFlow AML/MDS/MPD panel comprises two complementary marker combinations (tubes 29-32 and tubes 32-35; see Table 9), all of them containing HLA-DR-1, CD45-2, CD34-5 and CD117-6 as backbone markers (Table 9). The first set of tubes (tubes 29-32) is designed to preferably be used in combination with EuroFlow ALOT in order to aim at the detection and classification (lineage assignment and maturation) of myeloid malignancies, such as in AML and MDS, with a major focus on immature neutrophilic lineage (tube 30), monocytic lineage (tube 31), and erythroid lineage (tube 32). These four tubes also contribute to the detection of APL (tube 32), PNH (tubes 29 and 30) and other aberrant myeloid phenotypes.

25

The second set of tubes (tubes 33-35) provides additional information about megakaryocytic, basophilic, plasmacytoid dendritic lineages (tube 33), as well as

relevant information for the diagnosis of mastocytosis in association (or not) with AML/MDS (tube 34). Tube 35 further characterizes AML/MDS, and it is particularly focused on the aberrant expression of lymphoid-associated markers and abnormal lymphoid maturation.

5

TABLE 9. Multi-tube EuroFlow classification combinations for AML/MDS/MPD

Tube#	1	2	3	4	5	6	7	8	Aim
29	HLADR ^{BB}	CD45 ^{BB}	CD16	CD13	CD34 ^{BB}	CD117 ^{BB}	CD11b	CD10	Diagnosis of AML
30	HLADR ^{BB}	CD45 ^{BB}	CD35	CD64	CD34 ^{BB}	CD117 ^{BB}	IREM2	CD14	Diagnosis and subclassification of AML and PNH especially focussed on neutrophilic lineage
31	HLADR ^{BB}	CD45 ^{BB}	CD36	CD105	CD34 ^{BB}	CD117 ^{BB}	CD33	CD71	Diagnosis and subclassification of AML especially focussed on monocytic lineage
32	HLADR ^{BB}	CD45 ^{BB}	CD15	NG2	CD34 ^{BB}	CD117 ^{BB}	CD7	CD38	Diagnosis and subclassification of AML and APL especially focussed on erythroid lineage
33	HLADR ^{BB}	CD45 ^{BB}	CD42a plus CD61	CD203c	CD34 ^{BB}	CD117 ^{BB}	CD123	CD4	Diagnosis and subclassification of AML especially focussed on megakaryocytic, basophilic, and plasmacytoid dendritic lineages
34	HLADR ^{BB}	CD45 ^{BB}	CD41	CD25	CD34 ^{BB}	CD117 ^{BB}	CD42b	CD9	Characterization of AML/MDS
35	HLADR ^{BB}	CD45 ^{BB}	TdT	CD56	CD34 ^{BB}	CD117 ^{BB}	CD22	CD19	Characterization of AML/MDS

10 BB = Backbone markers; This multi-tube set of antibody combinations is not only useful for AML, MDS and MPD detection, but also for mastocytosis, PNH, ICUS. The tubes can also be successfully applied for disease staging and monitoring.

REFERENCES

- 1.- Consensual European Immunophenotyping panels for leukemia.
www.leukemia-net.org/content/home/.

- 5 2.- Wood BL, Arroz M, Barnett D, DiGiuseppe J, Greig B, Kussick SJ, Oldaker T,
Shenkin M, Stone E, Wallace P. 2006 Bethesda International Consensus
recommendations on the immunophenotypic analysis of hematolymphoid neoplasia
by flow cytometry: optimal reagents and reporting for the flow cytometric diagnosis
of hematopoietic neoplasia. *Cytometry B Clin Cytom.* 2007;72 Suppl 1:S14-22.
- 10

- 3.- Bene MC, Castoldi G, Knapp W, Ludwig WD, Matutes E, Orfao A, van't Veer
MB. Proposals for the immunological classification of acute leukemias. *European
Group for the Immunological Characterization of Leukemias (EGIL). Leukemia.*
1995;9:1783-6
- 15

- 4.- Rawstron AC, Orfao A, Beksac M, Bezdicikova L, Brooimans RA, Bumbea H,
Dalva K, Fuhler G, Gratama J, Hose D, Kovarova L, Lioznov M, Mateo G, Morilla
R, Mylin AK, Omedé P, Pellat-Deceunynck C, Perez Andres M, Petrucci M, Ruggeri
M, Rymkiewicz G, Schmitz A, Schreder M, Seynaeve C, Spacek M, de Tute RM,
20 Van Valckenborgh E, Weston-Bell N, Owen RG, San Miguel JF, Sonneveld P,
Johnsen HE; European Myeloma Network. Report of the European Myeloma
Network on multiparametric flow cytometry in multiple myeloma and related
disorders. *Haematologica.* 2008;93:431-8.

- 25 5.- Rawstron AC, Villamor N, Ritgen M, Böttcher S, Ghia P, Zehnder JL, Lozanski
G, Colomer D, Moreno C, Geuna M, Evans PA, Natkunam Y, Coutre SE, Avery ED,
Rassenti LZ, Kipps TJ, Caligaris-Cappio F, Kneba M, Byrd JC, Hallek MJ,
Montserrat E, Hillmen P. International standardized approach for flow cytometric
residual disease monitoring in chronic lymphocytic leukaemia. *Leukemia.*
30 2007;21:956-64.

- 6.- Matutes E, Owusu-Ankomah K, Morilla R, Garcia Marco J, Houlihan A, Que
TH, Catovsky D. The immunological profile of B-cell disorders and proposal of a
scoring system for the diagnosis of CLL. *Leukemia.* 1994;8:1640-5.

Claims

1. A reagent composition for flow cytometric immunophenotyping of leukocytes comprising at least eight distinct fluorochrome-conjugated antibodies comprising a set of at least three identification antibodies for the identification of a leukocyte population of interest and at least four characterization antibodies for further characterization and/or classification of said leukocyte population, wherein the antibodies are directed against the following combination of markers:

CD20, CD4, CD45, CD19, Ig λ , CD8, Ig κ , CD56, CD5, TCR $\gamma\delta$, CD3 and CD38,

wherein the antibody within each of the pairs CD20/CD4, Ig λ /CD8 and CD19/TCR $\gamma\delta$ is conjugated to the same fluorochrome, and wherein between different pairs the fluorochromes are distinguishable.

2. A set of at least two reagent compositions, said set comprising a reagent composition according to claim 1 and at least one further reagent composition comprising distinct fluorochrome-conjugated antibodies directed against one of the following combinations of markers:

- (i) CD20, CD45, CD23, CD10, CD79b, CD19, CD200 and CD43
- (ii) CD20, CD45, CD31, LAIR1, CD11c, CD19, IgM and CD81
- (iii) CD20, CD45, CD103, CD95, CD22, CD19, CXCR5 and CD49d
- (iv) CD20, CD45, CD62L, CD39, HLADR, CD19, CD27 and CD31
- (v) CD4, CD45, CD7, CD26, CD3, CD2, CD28 and CD8
- (vi) CD4, CD45, CD27, CCR7, CD3, CD45RO, CD45RA and CD8
- (vii) CD4, CD45, CD5, CD25, CD3, HLADR, cyTCL1 and CD8
- (viii) CD4, CD45, CD57, CD30, CD3, CD11c and CD8
- (ix) CD4, CD45, cyPerforin, cyGranzyme, CD3, CD16, CD94 and CD8
- (x) CD4, CD45, CD279, CD3 and CD8
- (xi) CD2, CD45, CD7, CD26, CD3, CD56, CD5 and CD19
- (xii) CD16, CD45, CD57, CD25, CD3, CD56, CD11c and CD19
- (xiii) HLADR, CD45, cyPerforin, cyGranzyme, CD3, CD56, CD94 and CD19; or
- (xiv) CD45, CD138, CD38, CD28, CD27, CD19, CD117 and CD81.

3. The set of reagent composition of claim 2, wherein each reagent composition comprises antibodies conjugated to pacific blue (PacB) or Horizon V450, pacific orange (PacO) or AMCA, fluorescein isothiocyanate (FITC) or Alexa488, phycoerythrin (PE), peridinin chlorophyl protein/cyanine 5.5 (PerCP-Cy5.5), PerCP or PE-TexasRed, phycoerythrin/cyanine7 (PE-Cy7), allophycocyanine (APC) or Alexa647, and allophycocyanine/H7 (APC-H7), APC-Cy7, Alexa680 or Alexa700.
4. The set of reagent compositions of claim 2 or 3, wherein each reagent composition comprises antibodies conjugated to pacific blue (PacB) or Horizon V450, pacific orange (PacO) or AMCA, fluorescein isothiocyanate (FITC) or Alexa488, phycoerythrin (PE), peridinin chlorophyl protein/cyanine 5.5 (PerCP-Cy5.5), PerCP or PE-TexasRed, phycoerythrin/cyanine7 (PE-Cy7), allophycocyanine (APC) or Alexa647, and allophycocyanine/H7 (APC-H7), APC-Cy7, Alexa680 or Alexa700.
5. The set of reagent compositions of claim 4, comprising at least a first reagent composition comprising an antibody against CD20 and an antibody against CD4 that are conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against Ig λ and an antibody against CD8 that are conjugated to FITC or Alexa488; an antibody against Igk and an antibody against CD56 that are conjugated to PE; an antibody against CD5 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD19 and an antibody against TCR $\gamma\delta$ that are conjugated to PE-Cy7; an antibody against CD3 that is conjugated to APC or Alexa647; and an antibody against CD38 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;
together with at least a second reagent composition selected from the group consisting of:
- (1) an antibody against CD20 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD23 that is conjugated to FITC or Alexa488; an antibody against CD10 that is conjugated to PE; an antibody against CD79b that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD19 that is conjugated to PE-Cy7;

an antibody against CD200 that is conjugated to APC or Alexa647; and an antibody against CD43 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(2) an antibody against CD20 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD31 that is conjugated to FITC or Alexa488; an antibody against LAIR1 that is conjugated to PE; an antibody against CD11c that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD19 that is conjugated to PE-Cy7; an antibody against IgM that is conjugated to APC or Alexa647; and an antibody against CD81 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(3) an antibody against CD20 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD103 that is conjugated to FITC or Alexa488; an antibody against CD95 that is conjugated to PE; an antibody against CD22 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD19 that is conjugated to PE-Cy7; an antibody against CXCR5 that is conjugated to APC or Alexa647; and an antibody against CD49d that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(4) an antibody against CD20 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD62L that is conjugated to FITC or Alexa488; an antibody against CD39 that is conjugated to PE; an antibody against HLADR that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD19 that is conjugated to PE-Cy7; an antibody against CD27 that is conjugated to APC or Alexa647; and an antibody against CD31 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(5) an antibody against CD4 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD7 that is conjugated to FITC or Alexa488; an antibody against CD26 that is conjugated to PE; an antibody against CD3 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD2 that is conjugated to PE-Cy7; an antibody against CD28 that is conjugated to APC or Alexa647; and an antibody against CD8 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(6) an antibody against CD4 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD27 that is conjugated to FITC or Alexa488; an antibody against CCR7 that is conjugated to PE; an antibody against CD3 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD45RO that is conjugated to PE-Cy7; an antibody against CD45RA that is conjugated to APC or Alexa647; and an antibody against CD8 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(7) an antibody against CD4 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD5 that is conjugated to FITC or Alexa488; an antibody against CD25 that is conjugated to PE; an antibody against CD3 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against HLADR that is conjugated to PE-Cy7; an antibody against cyTCL1 that is conjugated to APC or Alexa647; and an antibody against CD8 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(8) an antibody against CD4 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD57 that is conjugated to FITC or Alexa488; an antibody against CD30 that is conjugated to PE; an antibody against CD3 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD11c that is conjugated to APC or Alexa647; and an antibody against CD8 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(9) an antibody against CD4 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against cyPerforin that is conjugated to FITC or Alexa488; an antibody against cyGranzyme that is conjugated to PE; an antibody against CD3 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD16 that is conjugated to PE-Cy7; an antibody against CD94 that is conjugated to APC or Alexa647; and an antibody against CD8 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700 ;

(10) an antibody against CD4 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD279 that is conjugated to PE; an antibody against CD3 that is conjugated to

PerCP-Cy5.5, PerCP, or PE-Texas Red; and an antibody against CD8 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(11) an antibody against CD2 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD7 that is conjugated to FITC or Alexa488; an antibody against CD26 that is conjugated to PE; an antibody against CD3 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD56 that is conjugated to PE-Cy7; an antibody against CD5 that is conjugated to APC or Alexa647; and an antibody against CD19 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(12) an antibody against CD16 that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against CD57 that is conjugated to FITC or Alexa488; an antibody against CD25 that is conjugated to PE; an antibody against CD3 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD56 that is conjugated to PE-Cy7; an antibody against CD11c that is conjugated to APC or Alexa647; and an antibody against CD19 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(13) an antibody against HLADR that is conjugated to PacB or Horizon V450; an antibody against CD45 that is conjugated to PacO or AMCA; an antibody against cyPerforin that is conjugated to FITC or Alexa488; an antibody against cyGranzyme that is conjugated to PE; an antibody against CD3 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD56 that is conjugated to PE-Cy7; an antibody against CD94 that is conjugated to APC or Alexa647; and an antibody against CD19 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700;

(14) an antibody against CD45 that is conjugated to PacB or Horizon V450; an antibody against CD138 that is conjugated to PacO or AMCA; an antibody against CD38 that is conjugated to FITC or Alexa488; an antibody against CD28 that is conjugated to PE; an antibody against CD27 that is conjugated to PerCP-Cy5.5, PerCP, or PE-Texas Red; an antibody against CD19 that is conjugated to PE-Cy7; an antibody against CD117 that is conjugated to APC or Alexa647; and an antibody against CD81 that is conjugated to APC-H7, APC-Cy7, Alexa680, or Alexa700.

6. A diagnostic kit for flow cytometric immunophenotyping of leukocytes comprising a set of at least two reagent compositions according to any one of claims 2-5, optionally together with instructions for use, buffer, and/or control samples.
7. The diagnostic kit of claim 6 for the identification and characterization of mature lymphoid cells, comprising a set of reagent compositions recited in claim 2 under (i) through (v); (vi) through (x) ; (xi) through (xiii); or (xiv).
8. A method for flow cytometric immunophenotyping of leukocytes, comprising the steps of
 - (a) contacting a first aliquot of a biological sample comprising leukocytes with a first reagent composition of a set according to any one of claims 2 to 5 and contacting at least a second aliquot of said sample with a further reagent composition of said set;
 - (b) analyzing leukocytes in said aliquots in a flow cytometer; and
 - (c) storing and evaluating the data obtained.
9. The method of claim 8, wherein said sample is peripheral blood, bone marrow, tissue sample, or other type of body fluid.
10. The method of claim 9, wherein said sample is a cerebrospinal fluid, a vitreous fluid, a synovial fluid, pleural effusions or ascites.
11. The method of claim 9, wherein said tissue sample is lymph nodes, adenoid, spleen, or liver.
12. The method of claim 9, wherein said other type of body fluid is cerebrospinal fluid, vitreous fluid, synovial fluid, pleural effusions or ascites.
13. The method of any one of claims 8 to 12, wherein step (c) comprises the use of software for data integration and multidimensional analysis of flow cytometry files.

14. Method of claim 13, wherein said software is INFINICYT™.

Figure 1





Figure 2



EuroFlow panels of antibody combinations

Category 1

Mature lymphoid panel

Single tube screening combinations

- LST
- SST
- PCST

Multi-tube classification combinations

- B-CLPD
- T-CLPD
- NK-CLPD
- PCD

Category 2

Immature lymphoid panel

Single tube screening combination

- ALOT

Multi-tube classification combinations

- BCP-ALL
- T-ALL

Category 3

Myeloid panel

Single tube screening combination

- ALOT

Multi-tube classification combination

- AML/MDS/MPD*

* also suited for:

- mastocytosis
- PNH
- ICUS