

(12) PATENT
(19) AUSTRALIAN PATENT OFFICE

(11) Application No. AU 199538509 B2
(10) Patent No. 722176

(54) Title
Use of CYB medium for the transportation and storage of sperm

(51)⁷ International Patent Classification(s)
A01N 001/02

(21) Application No: **199538509**

(22) Application Date: **1995.11.09**

(87) WIPO No: **WO97/16965**

(43) Publication Date : **1997.05.29**

(43) Publication Journal Date : **1997.07.24**

(44) Accepted Journal Date : **2000.07.27**

(71) Applicant(s)
Medical Research Council; Applied Research Systems ARS Holding N.V.

(72) Inventor(s)
Robert John Aitken

(74) Agent/Attorney
F.B. RICE and CO.,139 Rathdowne Street,CARLTON VIC 3053

(56) Related Art
FERTILITY STERILITY, VOL. 35, NO.2, FEB 1981, PAGES 205-8
FERTILITY STERILITY, VOL. 34, NO. 6, DEC 1980, PP 607-9
JOURNAL OF ANDROLOGY, VOL.8, NO. 1, 1987, PAGES 41-7

OPI DATE 29/05/97 APPLN. ID 38509/95
AOJP DATE 24/07/97 PCT NUMBER PCT/GB95/02635



AU9538509

II

I)

(51) International Patent Classification ⁶ : A01N 1/02		A1	(11) International Publication Number: WO 97/16965
			(43) International Publication Date: 15 May 1997 (15.05.97)
(21) International Application Number: PCT/GB95/02635 (22) International Filing Date: 9 November 1995 (09.11.95) (71) Applicants (for all designated States except US): MEDICAL RESEARCH COUNCIL [GB/GB]; Centre for Reproductive Biology, Chalmers Street, Edinburgh EH3 9EW (GB). APPLIED RESEARCH SYSTEMS ARS HOLDING N.V. [NL/NL]; 14 John B Gorsiraweg, Curaçao (AN). (72) Inventor; and (75) Inventor/Applicant (for US only): AITKEN, Robert, John [GB/GB]; Medical Research Council, Centre for Reproductive Biology, 37 Chalmers Street, Edinburgh EH3 9EW (GB). (74) Agents: CHAPMAN, Paul, William et al.; Kilburn & Strode, 30 John Street, London WC1N 2DD (GB).		(81) Designated States: AU, CA, JP, US, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE). Published <i>With international search report.</i>	
(54) Title: USE OF CYB MEDIUM FOR THE TRANSPORTATION AND STORAGE OF SPERM			
(57) Abstract The use of CYB medium for the storage and/or transportation of semen at ambient temperature.			

USE OF CYB MEDIUM FOR THE TRANSPORTATION
AND STORAGE OF SPERM

5 The present invention relates to a novel method for storing and transporting sperm which ensures its viability.

10 Diagnosis of male infertility has become increasingly more complex over recent years in view of the increasing sophistication of andrological techniques. This has led to the creation of specialised andrology centres which possess the necessary expertise in techniques such as computerised image analysis, cell biology and biochemistry, and can thus offer the most up to date
15 diagnostic techniques.

However, a limitation on the growth of such specialised centres is that patients have to attend the centre in order to provide a semen sample immediately prior to
20 analysis. Fresh samples of semen are required since semen loses viability and functional competence when left in the presence of seminal plasma for any length of time, usually after 1 hour.

25 There thus exists a need to provide a method for storing/preserving semen samples at ambient temperature such that they can be sent to andrological centres for diagnostic testing, removing the need for patients to attend personally.

30 The present inventors have now found that semen can be stored at ambient temperature in a particular medium known as an egg yolk buffer.

Throughout this specification, unless the context requires otherwise, the word "comprise", or variations such as "comprises" or "comprising", will be understood to imply the inclusion of a stated element, integer or step, or group of elements, integers or steps, but not the exclusion of any other element, integer or step, or group of elements, integers or steps.

Disclosure of the Invention

In a first aspect, the present invention provides CYB medium when used in the storage and/or transportation of a semen at ambient temperature.

CYB medium is a known cryoprotectant medium first described by Wendel, L. and Prins, G.S., J. Androl. 8:41-47 (1987). It consists of the following components:

40% TES/TRIS buffer

30% sodium citrate/fructose solution

20% fresh egg yolk

1% pen-strep solution (10,000 IU)

TES = N' - Tris (hydroxymethyl)methyl-2-aminoethane sulphonic acid.

TRIS = Hydroxymethylaminomethane.

The TES/TRIS buffer contains 8.66 g TES and 2.06 g TRIS per 200 ml of water.


The sodium citrate/fructose solution contains 5.88 g sodium citrate and 4.0 g fructose per 200 ml of water.

The constituents of the CYB medium are mixed together and centrifuged at 600 g for 2 x 10 mins to remove particulate matter. The pH is adjusted to 7.4.

In a preferred embodiment, the medium also includes an antioxidant. Suitable antioxidants include α -tocopherol (vitamin E), catalase, glutathione and mannitol. It is believed that these components will act as free radical scavengers. A particularly preferred antioxidant is α -tocopherol, present at a concentration of 1mM.

In a second aspect, the invention provides a method for storing and/or transporting semen at ambient temperature which comprises the step of bringing the semen into contact with CYB medium. In this aspect of the invention, the semen is preferably diluted approximately 1:1 with the CYB medium.

In preferred embodiments of all of the above-noted aspects, the semen is human semen.



Preferred features of each aspect of the invention are as defined for each other aspect, *mutatis mutandis*.

The invention will now be described by way of the following examples which should not be construed as in any way limiting the invention.

5

EXAMPLE 1




(A) METHODS

10 (i) Sample Density Protocol

10 μ l of sperm suspension was added to 190 μ l of sperm diluting fluid (SDF; see below), mixed well, and both sides of the counting chamber were filled (improved

15

20



Neubauer ruling).

This was then left to settle.

- 5 The number of sperm in 5 large squares running diagonally across the slide was counted, each large square having 16 smaller squares.

- 10 When sperm heads lay on the border, only those that lay on the top and left side border were included.

The number of sperm in 5 large squares = million sperm per ml (10^6 /ml).

- 15 Sperm Diluting Fluid (SDF)

NaHCO₃ 50 g

Formalin 10 ml

Made up to 1 litre in distilled water.

20

(ii) **Scoring Motility**

- 25 10 μ l of sperm suspension is placed on a slide and covered with a 24x24mm coverslip. Using an eyepiece graticule with a square grid, an area is defined and motile sperm within that area are counted. Non-motile sperm within the area are counted. The field is moved and repeat process until at least 100 sperm have been counted. %motility (motile sperm as % of total sperm counted) is calculated.

30

(iii) **Percoll Prepared Sperm**

3 ml of 100% prepared percoll solution (see below) was placed in the bottom of a test tube. 3 ml of 50%

prepared percoll solution (see below) was carefully layered on top of the 100% percoll (50% percoll was layered 1ml at a time using a 1 ml pipette).

5 2 ml of semen sample was layered on top of the gradient column. This was then centrifuged at 1900 rpm for 20 mins (500 g).

10 Seminal plasma was removed from the top of the column (saved frozen).

15 The sperm was separated in bands (50% from the middle of the gradient at the 50/100 interface, and 100% from the bottom of the tube). This was resuspended in 7-10 ml of BWB (see below) and centrifuged at 1900 rpm for 5 mins.

20 The supernatant was removed and the sperm pellets resuspended in a known volume of BWB (500 µl- 1ml depending on recovery of sperm pellet).

Density was recorded and the sperm concentration was adjusted to $20 \times 10^6/\text{ml}$.

DISCONTINUOUS PERCOLL GRADIENTS

25

100% Percoll Solution - 100 ml

10 ml of 10 x Earle's Balanced Salts Solution
FLOW LABS, IRVINE, SCOTLAND

30

90 ml of percoll
PHARMACIA LKB BIOTECHNOLOGY AB, UPPSALA, SWEDEN

6 ml of albuminar 5%

ARMOUR PHARMACEUTICAL COMPANY, EASTBOURNE, ENGLAND

3 mg of sodium pyruvate

5 0.37 ml of sodium lactate

200 mg of sodium hydrogen carbonate (NaHCO_3)

50% Percoll Solution

10 100% percoll solution 1:1 diluted with BWB.

BWB PREPARATION

15 BWB Stock - Made up in 1 litre of distilled water.

	NaCl	5.54 g
	KCl	0.356 g
	CaCl_2 (dihydrate)	0.250 g
20	KH_2PO_4	0.162 g
	$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$	0.294 g

25 BWB - 200 ml

	NaHCO_3	420 mg
	Glucose	200 mg
30	Sodium pyruvate	6 mg
	Albuminar 5% (= 0.3% Final)	12 ml
	Sodium lactate	0.74 ml

35 Penicillin/streptomycin 2.0 ml
GIBCO

Hepes buffer 4.0 ml
FLOW LABORATORY, IRVINE, SCOTLAND

40 BWB stock 181 ml

45

(B) SAMPLE TREATMENT

5 (i) Patient semen samples were obtained in 30 ml sterile plastic sample containers. A 30 minute liquification period was allowed prior to recording sample volume, density and round cell count.

(ii) The samples were then split into two aliquots and treated as follows:

10

One aliquot was prepared on percoll gradients as detailed above.

15

The remaining aliquot was mixed with a fixed volume (4.0 ml) of CYB medium, thereby providing at least a 1:1 dilution for most semen samples. Sample motility was checked again before packaging in an insulated polystyrene box for despatch by a courier service over a period of 24 hours.

20

Aliquots prepared on Percoll gradients were analysed as follows:

(a) Acrosome Reaction Assay

25

Sperm samples were prepared according to the 50%/100% percoll gradient protocol, adjusting the sperm density to $20 \times 10^6/\text{ml}$.

30

200 μl of sperm was added to an equivalent volume of A23187 free acid to give a final concentration of 1.25 μM or 2.5 μM A23187.

35

The mixture was incubated at 37°C for 3 hours.

The mixture was washed at 500 g for 5 minutes and the supernatant removed and resuspended in BWB at $20 \times 10^6/\text{ml}$.

5 Motility of the sample was recorded.

50 μl of sperm suspension was added to 500 μl of Hypo-osmotic Swelling Medium and incubated for 1 hour at 37°C .

10

The suspension was centrifuged for 5 minutes at 500 g, the supernatant discarded and the pellet resuspended in 50 μl of ice cold methanol, i.e. final sperm concentration $20 \times 10^6/\text{ml}$.

15

10 μl was put onto each spot of a Hendley slide and was allowed to air dry.

20

This was then overlayed with Lectin-Fluorescein Isothiocyanate (2 mg/ml in PBS) and incubated for 15 minutes in the dark.

Excess lectin was then washed with PBS.

25

A drop of Citifluor was put on each spot of the Hendley slide and the slide was then covered with a coverslip for examination under a fluorescence microscope.

30

Hypo-osmotic Swelling Medium

7.35 g sodium citrate

13.51 g fructose

1 litre of distilled water.

(b) Sperm Penetration Assay

Sperm select 1:1 was diluted with BWB.

5 The diluted sperm select was then loaded into flat capillary tubes (200 μ m).

One end of the tube was capped with Critoseal (Mackay and Lynn, Edinburgh).

10

The open end of the sperm select capillary tubes was placed into the semen sample (50 μ l).

15

This was then incubated at 37°C for 30 min at an angle of 20° (approx).

20

The sperm select and Penetrak tubes were removed from the sperm, placed on a marked microscope tube and the number of sperm present at 1, 3 and 4.5 cm from the open end of the tube was counted (counting the number of sperm in 2 fields - using a x 40 objective).

Aliquots sent by courier

25

30

On arrival of the samples, the samples were analysed for motility loss before percoll preparation and diagnostic assays were repeated to determine if the transport diluent was capable of sustaining these important diagnostic parameters in the "field" situation. Experimental variation is reduced by the use of strict correlated protocols and by employing the same technician in the assessment of subjective assays such as the acrosome reaction test.

RESULTS

(a) Sperm motility

5 21 semen samples chosen from a random selection of donors were investigated. Motility before and after shipment in CYB medium was measured.

TABLE 1

10

15

20

25

30

35

Sperm Progressive Motility (WHO: a + b) in % motile on n = 21 donors	
Semen t = 0	semen/CYB t = 24
8	24
11	14
17	16
22	20
26	46
30	65
32	30
36	30
44	48
47	59
48	55
49	58
53	30
54	39
60	56
62	71
64	68
73	64
45	55
76	61
80	69
mean = 46 SD = 22	mean = 47 SD = 19

40

These results show that Sperm Progressive Motility at t = 24 h in CYB is on average well related to the response at t = 0.

(b) **Acrosome Reaction Assay**

For the same 21 samples of semen as noted above, the sperm percentage of viable cells was as follows:

TABLE 2

Ionophore A23817			
Semen t = 3		semen/CYB t = 24	
21		34	
25		55	
27		28	
37		19	
43		66	
57		42	
59		61	
60		56	
61		83	
64		62	
70		80	
74		79	
78		62	
78		82	
79		75	
79		82	
79		81	
81		87	
87		90	
mean = 61		mean = 64	
SD = 21		SD = 21	

These data show that the % cells viable, following treatment with Ionophore A23187 at t = 24 h in CYB, is, on average, well related to the response at t = 3.

(c) **Sperm Penetration Assay**

Sperm penetration (at t = 1.5 cm, 3 cm and 4.5 cm) at t = 0 and after t = 24 h in CYB was as follows:

TABLE 3

Sperm penetration on n = 20 donors					
1.5cm t=0	3cm t=0	4.5cm t=0	1.5cm t=24	3cm t=24	4.5cm t=24
251	34	2	166	27	8
174	19	1	24	4	3
121	6	0	35	7	0
224	43	8	250	47	10
110	19	2	85	18	1
260	54	9	255	110	8
196	20	3	358	74	6
275	45	8	86	26	3
270	91	25	284	99	28
180	19	3	300	117	15
2	0	0	3	1	0
6	2	0	5	21	0
28	2	0	34	5	1
192	45	9	186	43	7
92	18	1	101	21	2
114	11	0	109	12	1
210	21	2	215	23	3
300	40	7	270	38	8
34	3	0	39	5	1
170	37	5	156	32	6
mean = 160	mean = 26	mean = 4	mean = 36	mean = 36	mean = 6
SD = 93	SD = 23	SD = 6	SD = 111	SD = 36	SD = 7

There was no statistically significant difference between the values measured at t=0 and measured at t=24 hr in CYB.

Example 2

General Semen Analysis (other than motility) of Spermatozoa incubated in CYB for 24 hours

Semen liquefaction was assessed as being either non-existent, moderate or normal. In samples examined, the presence of mucous threads, ie incomplete liquefaction was also backed up by drawing the sample through a pipette to determine if the fluid flowed freely. After mixing 1:1 with CYB and incubation for 24 hours the assessment was repeated.

Sample Number	Semen t=0	Semen/CYB t=24hr
1	Normal	Normal
2	Normal	Normal
3	Moderate	Moderate
4	Normal	Normal
5	None	Moderate

Clearly, incubation in CYB has no significant effect on this observed property of sperm.

Abnormal form analysis was performed on sperm in semen and mixed 1:1 with CYB at t=0 hours and t=24 hours. The cells were fixed in standard formalin solution before analysis, the cells being observed using a x40 objective.

Sample Number	CYB t=0	CYB t=24	Semen t=0	Semen t=24	% diff t=0	% diff t=24
1	65	54	61	55	17	10
2	68	83	78	69	18	11
3	68	52	54	48	24	11
4	70	76	71	73	8	3
5	74	74	79	78	0	1

Clearly, incubation with CYB does not appear to alter sperm morphology. Any discrepancies can be assumed to be due to sampling error.

Measurement of Antisperm Antibodies in CYB

Samples were treated with the blood serum of patients previously known to have high antisperm antibody levels. Thus, antibody was transferred to the test samples, allowing a clear estimation of the effect of CYB on antibody detection after 24 hours.

The semen sample was allowed to liquefy for 30 minutes and was then prepared on mini percoll gradients by centrifugation at 600g for 5 minutes.

The semen was then washed twice in 0.4% BSA in Earles culture medium by centrifugation at 200g for 5 minutes.

The semen was then resuspended in 500 μ l of 0.4% BSA in Earles culture medium.

500 μ l of donor semen was incubated with 500 μ l of serum for 1 hour. After incubation, the sample was washed twice with 0.4% BSA in Earles culture medium and then resuspended with 1ml of culture medium and then mixed 1:1 with CYB.

For the MAR test the washed sperm were reconstituted into their seminal plasma before the MAR scoring. This involved mixing the semen with O+ blood preparation and IgG anti-sera. The IBBA test is amore sensitive and specific assay for the detection of anti-sperm antibodies.

	MAR t=0	MAR t=24 CYB	IBBA t=24 CYB
Control	Negative	Negative	IgA/IgG -ve
Stock IgG	100% positive	100% positive	IgA 50% +ve IgG 90% +ve
1 in 2 IgG	100% positive	100% positive	IgA 10% +ve IgG 90% +ve
1 in 10 IgG	100% positive	100% positive	IgA 10% +ve IgG 90% +ve

The IBBA results corresponded to the results gained by the serum samples in previous assays. Thus, the action of CYB does not appear to impair the detection of anti-sperm antibodies.

THE CLAIMS DEFINING THE INVENTION ARE AS FOLLOWS:

1. CYB medium when used in the storage and/or transportation of semen at ambient temperature.

5

2. CYB medium for use in the storage and/or transportation of semen at ambient temperature, the medium comprising 40% TES/TRIS buffer, 30 % sodium citrate/fructose solution, 20% fresh egg yolk, 1% pen-strep solution (10,000 IU).

10

3. A method for the storage and/or transportation of semen at ambient temperature which comprises the step of bringing the semen into contact with CYB medium.

15

4. A method as claimed in claim 3 wherein the semen is diluted 1:1 with CYB medium.

20

5. A medium according to claim 1 or claim 2, or the method as claimed in claim 3 or claim 4, wherein the CYB medium further comprises an antioxidant.

25


6. A medium or method as claimed in claim 5 wherein the antioxidant is α -tocopherol, catalase, glutathione or mannitol.

7. A medium or method as claimed in claim 6 wherein the antioxidant is α -tocopherol, present at a concentration of 1mM.

30

8. A medium or method as claimed in claim 6 or claim 7 wherein the semen is human semen.

9. A method for the storage and/or transportation of semen substantially as hereinbefore described with reference to any one of the foregoing Examples.



A 5x5 grid of dots forming the number 8. The dots are arranged in a pattern that resembles a digital display of the digit 8.

A 4x4 grid of dots forming the number 10. The dots are arranged as follows: Row 1: (1,1), (1,2), (1,3), (1,4); Row 2: (2,1), (2,2), (2,3), (2,4); Row 3: (3,1), (3,2), (3,3), (3,4); Row 4: (4,1), (4,2), (4,3), (4,4).

