



(86) Date de dépôt PCT/PCT Filing Date: 1995/03/15
(87) Date publication PCT/PCT Publication Date: 1995/09/21
(45) Date de délivrance/Issue Date: 2004/11/23
(85) Entrée phase nationale/National Entry: 1996/09/06
(86) N° demande PCT/PCT Application No.: CA 1995/000144
(87) N° publication PCT/PCT Publication No.: 1995/024936
(30) Priorité/Priority: 1994/03/15 (9405002.8) GB

(51) Cl.Int.⁶/Int.Cl.⁶ A61M 1/10, A61M 1/36, A61M 16/00


(72) Inventeurs/Inventors:
MUTCH, WILLIAM ALAN C., CA;
LEFEVRE, GERALD ROBIN, CA

(73) Propriétaire/Owner:
BIOVAR LIFE SUPPORT INC., CA

(74) Agent: SIM & MCBURNEY

(54) Titre : REGULATION DES SYSTEMES D'ASSISTANCE CARDIO-RESPIRATOIRE

(54) Title: CONTROL OF LIFE SUPPORT SYSTEMS



(57) Abrégé/Abstract:

The flow of a biological fluid to an organ is computer-controlled so that the natural variation of such flow is simulated. Specifically described are control of a blood pump flow output during CPB to mimic normal pulsatile blood flow from the heart and control of a ventilator output to mimic normal breathing of healthy lungs. A pattern of variation over time of instantaneous flow of a biological fluid to an organ of a mammalian species is established, a variable control parameter for regulation of flow of the biological fluid to the organ is generated in accordance with the pattern, and the flow of biological fluid to the organ is controlled in accordance with the variable control parameter.




INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁶ : A61M 1/10, 1/36, 16/00		(11) International Publication Number: WO 95/24936
A3		(43) International Publication Date: 21 September 1995 (21.09.95)
(21) International Application Number: PCT/CA95/00144 (22) International Filing Date: 15 March 1995 (15.03.95) (30) Priority Data: 9405002.8 15 March 1994 (15.03.94) GB (71) Applicant (for all designated States except US): UNIVERSITY OF MANITOBA [CA/CA]; 311 Administration Building, Winnipeg, Manitoba R3T 2N2 (CA). (72) Inventors; and (75) Inventors/Applicants (for US only): MUTCH, William, Alan, C. [CA/CA]; 89 Cordova Street, Winnipeg, Manitoba R3N 0Z9 (CA). LEFEVRE, Gerald, Robin [CA/CA]; 1118 Wolseley Avenue, Winnipeg, Manitoba R3P 0G9 (CA). (74) Agent: STEWART, Michael, I.; Sim & McBurney, Suite 701, 330 University Avenue, Toronto, Ontario M5G 1R7 (CA).		(81) Designated States: AM, AT, AU, BB, BG, BR, BY, CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GE, HU, JP, KE, KG, KP, KR, KZ, LK, LR, LT, LU, LV, MD, MG, MN, MW, MX, NL, NO, NZ, PL, PT, RO, RU, SD, SE, SI, SK, TJ, TT, UA, US, UZ, VN, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG), ARIPO patent (KE, MW, SD, SZ, UG). Published <i>With international search report. Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.</i> (88) Date of publication of the international search report: 7 December 1995 (07.12.95)

(54) Title: CONTROL OF LIFE SUPPORT SYSTEMS

(57) Abstract

The flow of a biological fluid to an organ is computer-controlled so that the natural variation of such flow is simulated. Specifically described are control of a blood pump flow output during CPB to mimic normal pulsatile blood flow from the heart and control of a ventilator output to mimic normal breathing of healthy lungs. A pattern of variation over time of instantaneous flow of a biological fluid to an organ of a mammalian species is established, a variable control parameter for regulation of flow of the biological fluid to the organ is generated in accordance with the pattern, and the flow of biological fluid to the organ is controlled in accordance with the variable control parameter.



2184992

TITLE OF INVENTION
CONTROL OF LIFE SUPPORT SYSTEMS

FIELD OF INVENTION

5 The present invention relates to medical life support systems, and, in particular, to the control of cardiopulmonary bypass pumps for open heart surgery and mechanical ventilators to lungs.

BACKGROUND TO THE INVENTION

10 During cardiopulmonary bypass surgery (CPB), the most common operation conducted in North America, the heart is stopped and the blood which normally returns to the right side of the heart passes through a pump and oxygenating system and is returned to the aorta, thereby
15 bypassing the heart and lungs. The flow of blood is essentially non-pulsatile with a low amplitude waveform having monotonous regularity.

 Although a common procedure (in excess of 400,000 open heart procedures per annum are conducted in North
20 America) and although tremendous strides have been made so that open heart surgery is safer for patients, the procedure is not without its dangers. Although the vast majority of patients have marked improvement in their cardiac functional status following their procedure, of
25 concern is the potential for damage to other organ systems which can occur due to the need for CPB.

 The following consequences have been identified with conventional non-pulsatile CPB, namely metabolic acidosis, interstitial fluid accumulation, elevated
30 systemic vascular resistance, arteriovenous shunting and impaired brain oxygenation. Of greatest concern is the potential for neurologic damage. Increasingly, well conducted prospective trials have demonstrated an alarming rate of post-operative neuropsychologic
35 disturbances following cardiac surgery. Recent studies have shown that up to 60 percent of patients undergoing open heart surgery have neuropsychologic deficits

SUBSTITUTE SHEET

following their operation, so that as many as 240,000 patients per annum may develop neurologic abnormalities following cardiac surgery. These disturbances are subtle but involve higher cognitive functions of the brain.

Mechanical ventilation of the lungs represents one of the major accomplishments of modern medicine and is one of the cornerstones upon which modern surgery and intensive care is based. Despite many major advances, mechanical ventilation is still associated with a number of alterations in respiratory function which causes morbidity and mortality in patients requiring this type of support. Inability to maintain gas exchange remains one of the major limiting factors with regard to life-support of critically ill patients. Even in healthy patients being ventilated during elective surgery, alterations in gas exchange can be demonstrated. These relate to collapse of small airways and alveoli. Prevention of these alterations would likely represent a major advance in management of all patients requiring ventilatory support. Conventional mechanical ventilation is monotonously regular in delivery of set tidal volume and respiratory rate.

The monotonous regularity of pumping of blood during CPB and of set tidal volume and respiratory rate of a mechanical ventilator is in contrast to the intrinsic spontaneously variable rhythms of heart rate, blood pressure and respiration, associated with a normal functioning heart as well as the considerable range of tidal volume and respiratory rate which a healthy individual demonstrates during breathing.

The article "Predictive Control by Physical Activity Rate of a Total Artificial Heart During Exercise" Maeda et al, ASAIO Transactions, vol. 34, No. 3, pages 480-484, describes a control system for a total artificial heart which will increase cardiac output in response to increased metabolic demands associated with exercise. An accelerometer which senses movement produces signals

AMENDED SHEET

which increase in proportion to increased activity associated with treadmill exercise. Measured changes in cardiac output in an exercised animal with its own heart are correlated with changes in the output of the accelerometer. Similar changes in accelerometer readings are then used to control the output of an implanted total artificial heart in an animal which is exercised to the same degree. Although the cardiac output changes depending on the intensity of the activity, there is no attempt to vary the beat to beat heart rate or beat to beat blood pressures at any given level of exercise.

SUMMARY OF THE INVENTION

In the present invention, the operation of a blood pump and mechanical ventilator are controlled to provide a flow of blood on the one hand and medical gases on the other which is varied in a manner that closely mimics the natural variation action of the heart and lungs and thereby overcomes some of the defects noted above. The invention is not applicable only to these two devices but is applicable to regulation in control of flow of any biological fluid to any organ. Although the existence of such variability in biological fluid flow is known, no one has heretofore taken such variability into account during the application of life support systems.

Blood normally is pumped in a monotonously regular non-pulsatile fashion or low amplitude pulsatile manner.

In one aspect of the present invention, a predetermined pattern of variations over time of instantaneous changes in flow of a biological fluid to an independently-functioning normal organ of a mammalian species first is generated. The mammalian species may be the human to whom the procedure is to be applied, another human or another mammalian species which is a model for a human, such as, a dog or a pig. The generated pattern may be an actual pattern determined from the mammalian species or may be a computer simulation of the known

AMENDED SHEET

2184992

4

variation in the flow. The generated pattern generally is provided with a sufficient number of determinations as to be representative of normal variation. Depending on the procedure involved, the pattern of variation may be established for the appropriate change in flow. For example, for control of blood pump during CPB, a pattern of variation over time of instantaneous blood pressure and heart rate of an independently-functioning healthy heart is established. For control of a ventilator device, a pattern of variation over time of instantaneous respiratory rate and tidal volume of independently-functional normal lungs is established.

A variable control parameter then is generated for regulation of flow of the biological fluid to an organ during controlled life support conditions in accordance with the pattern. This control parameter is most readily achieved by computer processing of the predetermined pattern of variation. In effecting such computer processing, the individual values of the parameters in the pattern and the peak-to-peak time interval between the individual values are recorded and analyzed. For example, for control of a CPB pump, each of the individual blood pressures for the pattern of instantaneous blood pressure and heart rate and the time interval (heart rate) between each of the individual blood pressures are recorded. For control of a ventilator, each of the individual respiratory rates and tidal volumes for the pattern of instantaneous respiratory rate and tidal volumes and the time interval between each of the individual respiratory rates and tidal volumes are recorded.

The variable control parameter generated in the procedure of the invention depends on the flow of biological fluid being regulated. In the case of the CPB pump, a signal is generated corresponding in value to an individually-determined blood pressure for a period of

AMENDED SHEET

2184992

5

time corresponding to the heart rate for the difference between the one individually-determined blood pressure and the next individually-determined blood pressure of the pattern. In the case of the ventilator, a signal is generated corresponding in value to an individually-determined respiratory rate and tidal volume.

In the present invention, the next step is to control the flow of biological fluid to the organ during controlled life support conditions in accordance with the variable control parameter. In this way, the flow of biological fluid to the organ is effected in accordance with the pre-established pattern of variation over time and hence mimics the natural flow of the biological fluid to the organ.

The manner of control of the flow of biological fluid to the organ depends on the biological fluid and the organ concerned. For example, in the case of the control of a CPB pump, a control voltage is generated corresponding in magnitude to the generated signal from the variation pattern and the control voltage is applied to the pump to provide an output of blood from the pump to the body during cardiopulmonary bypass of a pressure proportional to the magnitude of the signal for the period of time (peak-to-peak time interval). The steps of generating a signal, generating a control voltage and applying the control voltage to the pump then is repeated for each next individually-determined blood pressure of the pattern. Depending on the duration of the operation and the number of individual determinations in the pattern, it may be necessary to repeat these steps again for the pattern, reading either from the beginning or in the reverse direction. In this way, a pulsatile flow of blood from the pump is provided to the CPB patient which mimics normal pulsatile blood flow from a healthy heart.

Similarly, for the control of ventilating gas from a ventilator, a control voltage is generated

AMENDED S.F.P.

corresponding in magnitude to the generated signal from the variable pattern and the control voltage is applied to the ventilator device to provide an output of ventilating gas from the ventilator device of a respiratory rate proportional to the magnitude of the signal. The steps of generating a signal, generating a control voltage and applying the control voltage to the ventilation device are repeated for each next individually-determined respiratory rate of the pattern. In this way, a variable flow of ventilating gas from the ventilator device to the lungs of the body during controlled life support conditions is provided which mimics normal breathing of healthy lungs.

As noted above, the present invention is applicable not only to control of a CPB pump or a mechanical ventilator but also to any other operation or device involving this control of a biological fluid to any organ. For example, the principles of the invention may be used in intra aortic balloon counterpulsation (IABC), the technique used to support patients, usually following CPB, when they are unable to maintain adequate cardiac output, until enough heart function has returned to permit its discontinuation.

The principles of the invention may be employed to improve hemodialysis by introducing variability to the pumping to provide improved diffusion across the dialysis membrane by promoting better mixing of blood and avoidance of areas of relatively stagnant flow and thereby decreasing dialysis time.

In addition, the present invention may be employed with extracorporeal membrane oxygenation (ECMO), which is a modification of CPB in which bypass is instituted to support the patient while giving the lungs a chance to heal. The patient is ventilated while on ECMO and if the therapy is successful, eventually weaned off ECMO and the ventilator. Computer control of the CPB pump in this

ALLEN 07/11/97

6A

situation has the potential to enhance organ perfusion, while computer controlled ventilation has the potential to activate lung healing and of improving gas exchange in order to facilitate earlier weaning from ECMO.

The present invention further may be employed in conjunction with right and left ventricular assist devices (RVAD and LVAD), which are occasionally used to support patients after CPB when they are unable to maintain adequate output without this type of support. Patients given such support are simultaneously being ventilated. Computer-controlled ventilation and computer-controlled RVAD and LVAD, as provided herein, may improve organ perfusion while computer-controlled ventilation may indirectly influence hemodynamic variability.

Another application of the principles of the present invention is in the perfusion of organs prior to transplantation.

In accordance with another aspect of the present invention, there is provided apparatus for controlling the flow of a biological fluid to an organ, which comprises: means for establishing a predetermined pattern of variations over time of instantaneous changes in flow of a biological fluid to an independently-functioning normal organ of a mammalian species; means for generating a variable control parameter for regulation of flow of the biological fluid to an organ during controlled life support conditions in accordance with the predetermined pattern and means for controlling the flow of the biological fluid to the organ during controlled life support conditions in accordance with the variable control parameter.

When the apparatus provided in accordance with this aspect of the present invention is used in controlling the flow of blood by a pump to a body during

PENDING SHEET

cardiopulmonary by-pass wherein said means for establishing a predetermined pattern comprises means for establishing a predetermined pattern of variation over time of instantaneous blood pressure and heart rate of an independently-functioning heart of a mammalian species; said means for generating a variable control parameter comprises means for generating a signal corresponding in value to an individually-determined blood pressure for a period of time corresponding to the heart rate for a difference between such individually-determined blood pressure and the next individually-determined blood pressure of said predetermined pattern and means for generating a control voltage corresponding in magnitude to said signal; and said means for controlling the flow of the biological fluid comprises means for applying said control voltage to said pump to provide an output of blood from said pump to the body during cardiopulmonary bypass of a pressure proportional to the magnitude of the signal for said period of time and means for repeating the steps of generating a signal, generating a control voltage and applying the control voltage to the pump for each next individually-determined blood pressure of said predetermined pattern.

When the apparatus provided in accordance with this aspect of the present invention is used in controlling the flow of ventilating gas from a ventilator device to the lungs of a body during controlled life support conditions wherein said means for establishing a predetermined pattern comprises means for establishing a predetermined pattern of variation over time of instantaneous respiratory rate and tidal volume of the independently-functioning healthy lungs of a mammalian species; said means for generating a variable control parameter comprises means for generating a signal corresponding in value to an individually-determined

6C

respiratory rate and tidal volume in said predetermined pattern and means for generating control voltage corresponding in magnitude to said signal; and said means for controlling the flow of biological fluid comprises means for applying the control voltage to the ventilator device to provide an output of ventilating gas from the ventilating device of a respiratory rate proportional to the said signal and means for repeating the steps of generating a signal, generating a control voltage and applying the control voltage to the ventilator device for each next individually-determined respiratory rate of the predetermined pattern.

BRIEF DESCRIPTION OF DRAWINGS

Figure 1 is a flow diagram showing the various steps of the operation of a CPB pump in accordance with one embodiment of the invention.

RECEIVED SHEET

Figure 2 shows a typical plot of natural variation of systolic blood pressure over time. In this instance, in a dog, data is captured to a data acquisition system from a dog. Following processing, these data are used as
5 an input file for the computer controller used to vary roller pump head revolutions/min based on variability in beat-to-beat intervals and pressure.

Figure 3 shows a typical pump pulsation profile for a blood pump controlled in accordance with the invention
10 superimposed upon a plot of peak-to-peak variations in blood pressure derived from a plot similar to that shown in Figure 2.

Figure 4 shows typical circuitry for computer control of a blood pump motor according to the invention.

15 Figure 5 shows a wiring diagram for a blood pump motor.

Figure 6 shows the changes in respiratory rate (breaths/min) over time. Such data is used to create the input file for the ventilator computer controller. These
20 data have a mean rate of 15 breaths/min.

Figure 7 shows the arterial minus superior sagittal sinus (cerebral venous) oxygen content difference during CPB with computer-controlled and conventional roller pump bypass in dogs (n = 6 both groups). The oxygen content
25 difference is stable during conduct of CPB in the computer-controlled group (operation in accordance with the invention). In the control group (conventional roller pump), increased oxygen extraction occurs during the period of rewarming compared to the hypothermic
30 period of CPB immediately before. Such changes in oxygen content difference with rewarming are associated with cognitive impairment following CPB in man.

Figure 8 shows a comparison of arterial oxygenation (PaO_2 in mm Hg) for computer-controlled vs. conventional
35 mechanical ventilation in a porcine oleic acid lung injury model. Significantly greater PaO_2 is seen at the

SUBSTITUTE SHEET

latter time periods in the experiment in animals ventilated with the computer-controlled ventilator (operation in accordance with the invention).

Figure 9 is similar to Figure 6 and shows the changes in respiratory rate (breathes/min) over time as well as the mean value.

Figure 10 shows a frequency vs respiratory rate plot devised from the graph of Figure 9.

Figure 11 is a block diagram showing the overall signal flow to the ventilator (Ohio 7000) used in the Examples described below. The 'Ohio Interface' module is connected to the Metrabyte DASH16 analog to digital (A/D) converter via bi-directional control lines via the 'DAS16 Cable' and 'DAS16 Jumper Box'.

Figure 12 is the internal wiring harness of the Ohio Interface Unit. Connector 'P1' and module 'Ohio 7000' refer to the cable and modifications added to the Ohio 7000 Ventilator. This Figure shows all physical wiring connections of the electronic modules ('Volume Modulator' and 'Rate Modulator') to the switches and connectors. The Ohio 7000 supplies the power (8V & COM). Connectors 'P2' and 'P3' are opposite gender. Experimental monitoring jacks ('J1' through 'J4') are for an external data acquisition system. Connector 'P3' interfaces to the Metrabyte model DASH16 A/D and digital to analog (D/A) converter. Switch S1 and S2 provide the ability to cancel 'RATE' or 'VOLUME' modulation individually ('manual' position) at any time.

Figure 13 displays an excerpt of the Ohio 7000 Ventilator control circuitry generated with an ohm meter and a photocopy of the service manual. Refer to the 'IE' ratio control R14. The control's wiper was originally connected to amplifier U22B. Modulation of the 'IE' control is introduced by inserting an external summing amplifier between the control and amplifier U22B (pins 3 and 7 of connector P4). Similarly, the 'RATE' and

SUBSTITUTE SHEET

'VOLUME' controls and their associated amplifiers (U22A and U5D) are routed to P4 pins 4 and 8, and 5 and 9 respectively. The 2 volt reference voltage for the 'RATE' and 'VOLUME' controls, as well as 8 volt supply and common are also routed to connector P4. Referring back to connector P2 in Figure 12: external modulation of 'IE' was disabled (pins 3 and 7 jumpered), and the 2 volt reference (pin 2) was not required.

Figures 14 and 15 display the modulation control electronics for 'VOLUME' and 'RATE' respectively. The circuit boards are interchangeable. Referring to Figure 14, the wiper of the 'VOLUME' control is buffered by amplifier U1D, and coupled to a non-inverting summing amplifier (U1C), whose output is routed back to the 'VOLUME' control's original destination (label 'VOL CTRL' is connected to Figure 13 amplifier U5D). The 'VOLUME' output at U1D is also routed to the A/D converter channel 0 input (U2A), and the 'VOL POT MONITOR' (U2B), which allows optional monitoring of the Ohio 7000's 'VOLUME' control level before modulation is inserted (J1 on Figure 2). The modulation reference signal from the D/A converter channel 0 output is buffered and scaled (I1A and U1B) and routed to the other input of the non-inverting summing amplifier (U1C), whose output is also routed to A/D converter channel 1 input (U2C), and to U2D - the 'VOL MONITOR' (J2 on Figure 12), which allow optional monitoring of 'VOLUME' after modulation is inserted. Since the D/A converter's output is in a range between 0 and 5 volts, the modulation level only increases the 'VOLUME', hence the position of the ventilator's 'VOLUME' control sets the minimum or baseline level of 'VOLUME'.

Figures 16 and 17 document (Figure 6) the cable used to interface the Ohio Modulation Unit to the 'DAS16 Jumper Box' (Figure 7), which is connected to the 37 pin connector on the Metrabyte DASH16 A/D and D/A converter.

SUBSTITUTE SHEET

DESCRIPTION OF PREFERRED EMBODIMENT

In the following description of a preferred embodiment, there is description of the application of the present invention to control of a blood pump.

5 However, it will be understood that the principles described with respect to such blood pump embodiment apply to other devices, including control of ventilators as described elsewhere. The steps involved are shown schematically in Figure 1.

10 During a CPB procedure, an electrically-driven pump is used to maintain blood flow, as described above. Generally, a roller pump is employed for this function, in which a pair of diametrically opposed rotating arms engage a flexible tube through which blood is forced by
15 the action of the arms engaging the flexible tube.

An input file for a computer-controller for the pump first is established for the variation of systolic pressure with time for a typical animal, such as a human, a dog or a pig. A typical plot of the gross variation of
20 systolic pressure in mm Hg over time is shown in Figure 2.

This information, which may contain many thousands of observations of systolic pressure, is loaded into the computer memory and processed by peak height analysis.
25 In this peak height analysis, the maximum, minimum average systolic pressures are determined and may be displayed, the minimum values are removed and the minimum, maximum and average of the remaining peak pressures is recalculated and, if desired, displayed.
30 This information then determines the pulse pressure amplitude and beat-to-beat heart beat variation in the pattern.

The lowest and highest values of the pressure values from the peak-to-peak variation then are established and
35 these values are used to set a minimum and maximum rate for the blood pump, respectively, which then determines

SUBSTITUTE SHEET

the maximum amount of computer modulation. For example, a baseline pressure of 80 mm Hg with a 20 mm Hg variation may be established based on the input file, which then provides a peak pressure ranging from 80 to 100 mm Hg.

5 The computer digital output signal corresponding in magnitude to a peak-to-peak value above the minimum is connected through a digital-to-analog (D/A) converter, which produces an analog voltage control signal to the blood pump to increase the blood pump rate. The computer

10 generates a voltage on the D/A converter proportional to the peak pressure variations for a time proportional to the beat-to-beat interval. The voltage then is used to increase the rpm of the pump from the minimum or baseline setting. The data stored in memory is converted into

15 time steps and relative amplitudes from 0 to 100%. For each time step, the D/A drive is held at the relative level until the next time step occurs. A typical pump pulsation profile superimposed upon a plot of peak-to-peak variations in blood pressure is shown in Figure 3.

20 As may be seen in Figure 3, over a period of approximately 18 seconds, the computer-controller output varies between 1 and 3.5 volts. The changes in roller pump speed has resulted in escalations in blood pressure varying between 82 and 102 mm Hg. The data stored in

25 memory is initially scanned in a forward direction for observations 1 to N. As necessary, the data is reverse scanned continuously from observations N-1 to 1 and then forwards from 2 to N. etc. until the program is terminated.

30 Any convenient form of the pump controller may be employed to receive the control signals from the computer and the corresponding voltage signals from the D/A converter. A typical circuitry is shown in Figure 4 while a wiring diagram appears as Figure 5. In this

35 circuit, a non-inverting summing amplifier with input buffering is provided, power to operate the amplifier is

SUBSTITUTE SHEET

2184992

12

from the roller pump rate controller. The signal from the original speed control is buffered by a buffer (IC-1), which is applied to one input of a summing amplifier (IC-2). The other input of the summing amplifier is
 5 received from the computer modulation signal received from the D/A converter via an external scaling box through buffer IC-3 and amplifier IC-4, which permits an increased voltage range, according to the desired multiple of the amplification, to be applied to the pump
 10 servo motor than provided by the D/A converter.

In this way, the roller pump revolutions/min are altered to recreate the pattern of spontaneous biologic variability in the heart function.

The computer operations described above may be
 15 effected on any convenient computer hardware programmed in any desired manner to effect the analysis described above to provide the blood pump control voltage. A program which may be employed, named Purfus, has the listing appearing in Table 1 below.

20 A configuration file, named Purfus Cfg, is necessary for the operation of Purfus program. This file contains a number corresponding to the base address of the D/A converter card:

	Decimal	Hex
25	768	300 H
	784	310 H
	800	320 H

Computer control of a ventilator to apply the principles of the invention thereto may be effected using
 30 the computer program shown in Table 2 below. This software allows the 'RATE' and the 'VOLUME' settings of the ventilator controls to be modulated independently via a data file, generated before hand and in a form such as appears in Figure 6.

35 In order to implement computer control of the ventilator using the circuitry shown in Figures 11 to 17,

SUBSTITUTE SHEET

a means of converting voltage to ventilation 'RATE' and 'VOLUME' is provided. A linear regression analysis of the ventilator's 'RATE' and 'VOLUME' potentiometers output voltage versus dial calibrations is performed ($R^2 = 0.9996$). Functions are converted 'RATE' and 'VOLUME' into voltage, and vice-versa. By control loop scanning the A/D converter channel in the background using the high speed DMA facility, the current setting of the ventilator's 'RATE' control is acquired. If the modulation level for the current time step is greater than the baseline (set by the current setting of the ventilator's 'RATE' control), the D/A converter channel generates a voltage level, which is passed to the summing amplifier in the Ventilator Modulation Unit, necessary to increase the current baseline value of 'RATE' to the modulation level of 'RATE'. The output of the summing amplifier is then sampled by another A/D converter channel, converted into 'RATE', and displayed on the computer screen. The 'VOLUME' are updated in each 'loop' of the control program which executes every 400 milliseconds on a 'control' computer (a 4.77 MHz 8088 processor with a 8087).

EXAMPLES

Example 1:

This Example illustrates the methods and materials used in ventilation experiments carried out on pigs.

Pig preparation:

Thirteen (13) pigs weighting 20 to 30 kg were studied. All pigs received atropine 0.6 mg and ketamine 10 mg/kg intramuscularly for induction of anesthesia. Once sedated, isoflurane in oxygen was administered by face mask. When airway reflexes had been obtunded, the pig was intubated with a 6.0 mm endotracheal tube. Mechanical ventilation was instituted with an Ohio 7000 anesthesia ventilator at 15 breaths/min with the minute ventilation adjusted to maintain the end-tidal CO_2 at 35

SUBSTITUTE SHEET

to 40 mm Hg. Isoflurane was administered at 2.0 percent end-tidal in oxygen during surgical preparation. Lactated Ringer's was infused IV at 10 ml/kg/hr during the experiment. Pancuronium bromide was administered IV
5 intermittently for muscle relaxation.

The animal was turned supine and a cutdown performed in the groin. A double-lumen catheter was placed in the femoral artery for intermittent sampling of blood for arterial blood gases (ABG) and continuous recording of
10 arterial pressure. A 7.5 Fr pulmonary artery catheter was inserted via the femoral vein and advanced with the balloon inflated until a pulmonary capillary wedge pressure (PCWP) was obtained. Pulmonary artery pressure was continuously recorded. Mixed venous blood was
15 sampled from the distal end of the pulmonary artery catheter. Cardiac output (CO) was measured intermittently, by thermodilution, following 5 ml injection of room temperature saline (performed in triplicate). Following surgery, the animal was allowed
20 to stabilize for 30 minutes and the isoflurane concentration was reduced to 1.5 percent end-tidal.

Baseline hemodynamic and respiratory measurements were then obtained. These included measurements of mean arterial pressure (MAP), mean pulmonary artery pressure
25 (MPAP), PCWP, airway pressures at the proximal end of the endotracheal tube (all recorded to a Gould 2600 oscillograph and to an advanced CODAS data acquisition system), and CO. Gas measurements included arterial and mixed venous blood gases and end-expired gas sampled from
30 the expiratory limb of the anesthesia circuit. These were measured using a Radiometer ABL3. Arterial and mixed venous oxygen content, oxygen saturation and hemoglobin concentration were measured with a Radiometer OSM3 set for porcine blood. All measurements were
35 obtained in duplicate. Calculated indices included

SUBSTITUTE SHEET

15 2184992

pulmonary vascular resistance (PVR), dead space ventilation (VD/VT) and shunt fraction (QS/QT).

Oleic Acid Lung Injury: After the above measurements were obtained a Valsalva maneuver was done (mean airway pressure 30 cm H₂O for 5 seconds). An infusion of oleic acid was started at 0.2 ml/kg/hr through the infusion port of the pulmonary artery catheter. At 5 min intervals the Valsalva maneuver was repeated and 1 min later an arterial blood gas obtained. The oleic acid infusion was continued until the PaO₂ decreased to $\leq$ 200 mm Hg for 2 consecutive measurements. At this point the infusion was stopped and the volume infused noted. Following repeat hemodynamic and respiratory measurements as above, the animals were randomly allocated to one of two ventilatory modes; conventional IPPV with the respiratory rate (RR) fixed at 15 breaths/min with the minute ventilation (MV) changed to maintain PaCO₂ at $\leq$ 45 mm Hg (control), or IPPV with a computer-controller with variable RR but with a mean of 15 breaths/min (computer). Again, MV was adjusted to maintain PaCO₂ at $\leq$ 45 mm Hg. Ventilation continued with either the control or computer mode for the duration of the experiment. Every 30 minutes for 180 min, hemodynamic and respiratory data was obtained as above in duplicate. At 180 min, airway pressure data was acquired to the data acquisition system over a 2 min time period to sample approximately 30 consecutive breaths.

Computer-controlled Ventilation: The ventilation used in this experiment was an Ohio 7000 having electronics as illustrated in Figures 11 to 17 described above. The computer-controller software (Table 8) allows the rate and volume settings of the ventilator controls to be modulated independently via a previously generated data file. Data from this modulation file is depicted graphically in Figure 9. A frequency vs respiratory rate plot is shown in Figure 10. The modulation file was

SUBSTITUTE SHEET

2184992
16

generated from hemodynamic and respiratory excursions from an anesthetized dog. This information was captured by data acquisition, processed and scaled to produce breath-to-breath variability. Ventilatory variability
5 can also be measured directly, stored and subsequently used to control the ventilation. Hardware was also developed to implement computer-control of the ventilator which necessitated converting voltage outputs from the breath-to-breath variability period to control ventilator
10 respiratory rate and tidal volume, as seen in Figures 11 to 17. In this experimental configuration, only RR was changed. As there was employed a ventilator which functioned as a volume divider, change in the RR resulted in reciprocal changes in the TV. Functions were
15 developed to convert ventilator rate and volume into voltage and vice versa. Output to control RR was updated every 400 msec and changed accordingly based on the modulation data file. The computer ventilator RR was set to 10 breaths/min baseline and the modulation file
20 programmed ventilation from 10 to 22 breaths/min with a mean value of 15 breaths/min.

Post-hoc analysis: The data file of airway pressures was processed to integrate the area under the pressure time curve to give mean airway pressure. Mean peak
25 airway pressure was also calculated. Because of the variability in RR and TV in the computer-controlled ventilator group, a minimum of 25 breaths were analyzed in each experiment. At the end of each experiment, the animal was killed with a lethal dose of thiopental, and
30 a sternotomy done to remove the lungs. The lungs were weighed wet and then suspended and aerated to commence drying. The following day, the lungs were placed in an oven to dry to a stable weight ($\pm$ 5 percent on consecutive days). The wet:dry lung weight ratio was
35 calculated.

SUBSTITUTE SHEET

17 2184992

Statistical Analysis: Multiple comparisons of data within and between groups was with repeated measures ANOVA. A p-value ≤ 0.05 was considered significant for group x time interactions or differences between groups.

5 Least squares means test matrices were generated for post-hoc comparisons. Bonferroni's correction was applied when multiple comparisons were examined within groups. Single comparisons between groups were by Student's t-test, $p \leq 0.05$ considered significant.

10 Example 2

This Example provides the results of the experimentation described in Example 1.

The computer-controlled ventilator varied respiration from 10 to 22 breaths/min (mean $\pm$ SD; $15.0 \pm$
15 2.3). There were 369 RR and TV combinations over 1089 sec before the modulation file looped to repeat itself.

The demographic data from the experiments is shown in Table 3 below. There were 7 animals in the computer group and 6 in the control ventilator group. The animals
20 in the two groups did not differ for weight or in the amount of oleic acid infused to induce the lung injury. The mean airway pressure did not differ between groups nor did the mean peak airway pressure. There was no difference between groups in the wet:dry lung weight
25 ratio.

There was no difference between groups for blood or nasopharyngeal temperature (group x time interaction; $p = 0.1772$ and 0.2665 respectively) (Table 4 below). A group effect was seen for baseline blood temperature of
30 0.6 degrees. In both groups, temperature increased significantly following lung injury. A marked difference was seen between groups for hemoglobin concentration ($p = 0.0014$ group x time interaction). In both groups hemoglobin increased significantly following lung injury,
35 but continued to increase in the control group. There was no interaction for pH between groups ($p = 0.2325$) but

SUBSTITUTE SHEET

there was a group effect with lower pH in the latter periods of the experiment in the control group.

Hemodynamic data is shown in Table 5 below. The MAP was stable between groups (group x time interaction; $p = 0.4429$). In both groups MAP decreased significantly following lung injury. The MPAP showed an interaction ($p = 0.0198$). In both groups the MPAP increased markedly following oleic acid. Baseline MPAP was significantly higher in the computer group then lower by 90 minutes. The PCWP was essentially identical between groups. No interaction was seen for PVR but a marked group effect was seen ($p = 0.0001$). In both groups PVR increased dramatically with lung injury. The PVR was significantly higher from 90 min on in the control group. There was no difference between groups for cardiac output at any time period. In both, the CO decreased to about 60 percent of control and remained unchanged.

Respiratory gas data is shown in Table 6 below. End-expired CO_2 (PeCO_2) did not differ between groups. There was a significant increase in PaCO_2 following lung injury in both groups. This correlated to the significant increase in dead space ventilation (VD/VT) seen. Importantly, a significant group x time interaction was seen for PaO_2 ($p = 0.0448$). A markedly significant group effect was seen as well ($p = 0.0001$). This is evident from significantly greater PaO_2 at time periods 60 - 150 min after oleic acid infusion. Of note, at baseline and at Time 0, PaO_2 values are not significantly different. In both groups the shunt fraction (QS/QT) increased significantly following lung injury.

Example 3

This Example discusses the results obtained in the experiments described in the preceding Examples.

In these experiments described in Example 2, it has been demonstrated that oxygenation is improved by

SUBSTITUTE SHEET

modifying mechanical ventilation to incorporate biologic variability. Through use of a computer-controller, variability in RR and TV resulted in significantly improved PaO₂ compared to standard IPPV with the same ventilator. This improvement in oxygenation was accomplished without an increase in mean airway or mean peak airway pressures.

No differences were seen between the two groups for amount of oleic acid administered to injure the lungs. The wet:dry lung weight ratios suggest a similar injury between the two groups. The two groups were very similar at baseline and following lung injury for PCWP, CO and PVR. Similar increases in shunt fraction and dead space ventilation were also seen for these two groups over the same time periods. Thus, the two groups appear not to differ prior to being randomized to control or computer-controlled ventilation.

As configured for this study, the ventilator delivered 369 different RR and TV combinations with a mean RR of 15 breaths/min over 1089 sec (18.2 min). This is contrasted to a single RR of 15 breaths/min throughout the duration of the experiment in the control group. Some variability exists in the control group because MV was adjusted to attempt to maintain PaCO₂ ≤ 45 mm Hg when VD/VT increased with lung injury. However, this entailed only a maximum of 6 changes in MV, over the course of any one experiment, when MV was changed if PaCO₂ was out-of-range, at the end of any 30 min measurement period. Thus, over a 30 min measurement period, RR and TV were essentially fixed in the control group but there were 369 x 30/18.2 = 608 different RR and TV combinations in the computer group. If, as Suki et al. (Nature, Vol. 368, April 1994, p. 615-618) suggest, airway recruitment is stochastic, then the probability of airway opening is dramatically improved using the computer-controlled ventilator. The experimental results provided herein

SUBSTITUTE SHEET

20 2184992

indirectly suggest that this finding is so. Of greatest significance is that mean PaO_2 improved following lung injury in the computer group out to 150 min (Figure 8) compared to an inexorable decline in PaO_2 in the control group. This improved PaO_2 was associated with significantly lower MPAP and PVR at identical PCWP in the computer group suggesting better ventilation/perfusion matching and lower pulmonary vascular resistance at similar cardiac outputs. The marked difference in hemoglobin concentration following lung injury is an independent marker that the two ventilatory modes differ. The increasing hemoglobin concentration in the control group suggests further accumulation of lung water. In the computer group, hemoglobin concentration remained essentially stable. This is especially so as the initial mean increase in hemoglobin concentration is identical in the two groups (26 percent). At Time 0, both groups were control mode ventilated. Only after Time 0, was computer ventilation initiated. Therefore, by inference, during the conduct of the experiment, lung water accumulation was less in the computer group with associated better oxygenation.

Example 4

This Example illustrates the materials and methods used to evaluate a CPB pump in dogs, using the computer control operation described above with reference to Figures 1 to 4.

Experimental Preparation: Twelve mongrel dogs (21 ± 3 kg) were studied. All animals were anesthetized with sodium thiopental ($25 \text{ mg} \cdot \text{kg}^{-1}$). The trachea was intubated and the animal ventilated with O_2 . The minute ventilation was adjusted to maintain PaCO_2 at 35 to 40 mm Hg. The dog was positioned in a stereotactic head-frame in a modified sphinx position. Bipolar EEG electrodes were placed over the parietal hemisphere bilaterally and monitored by an Interspec Neurotrac® in raw EEG mode.

SUBSTITUTE SHEET

Temperature was measured in the nasopharynx using a calibrated YSI tele-thermometer®. Anesthesia was maintained with isoflurane 1.3% end-tidal (1 MAC) during the surgical preparation. Following thoracotomy, the isoflurane was discontinued for a minimum of 30 min and the EEG made isoelectric with a bolus of thiopental. A continuous infusion of thiopental was initiated at 10 mg.kg⁻¹.hr⁻¹ to maintain the EEG isoelectric during CPB. Neuromuscular relaxation was achieved with pancuronium bromide.

A flow-directed catheter was advanced through the left femoral vein into the right atrium for central venous pressure (CVP) monitoring. A femoral artery catheter was advanced into the distal aorta for arterial pressure (MAP) monitoring. A double lumen (7.5 FR) catheter was inserted into the left brachial artery for intermittent blood withdrawal. The superior sagittal sinus (SSS) was exposed by trephine and the posterior one-third was cannulated non-occlusively by insertion of a 22-gauge intravenous catheter. Continuous cerebrospinal fluid pressure (CSFP) measurements were recorded by inserting a 22-gauge spinal needle into the cisterna magna with the use of a micromanipulator (Narishige®). A right thoracotomy was performed. The right atrium and proximal aorta were cannulated with a single stage 38 Fr atrial and Jostra® 21 Fr or 24 Fr aortic cannula, respectively. Following the initiation of CPB, the left ventricle was vented by a cannula inserted through the right superior pulmonary vein and the proximal aorta was occluded with a Seldinger vascular clamp.

All blood pressures and the CSFP were measured by calibrated Abbott® transducers referenced to the intra-auricular line. Data were recorded continuously on paper by an oscillograph (recorder model 7754A®, Hewlett Packard) and intermittently on hard disk by an IBM PC-AT®

SUBSTITUTE SHEET

computer based data acquisition system (Dataq Instruments®). The latter data are reported. Arterial and SSS blood gases were measured before and after each microsphere injection by an ABL-3 Acid-Base Laboratory (Radiometer®) at 37°C and not corrected for temperature. Arterial and cerebral venous (SSS) oxygen content and hemoglobin were measured by Radiometer OSM-3 (specific for canine blood).

Cardiopulmonary bypass was conducted utilizing a Travenol® non-pulsatile roller pump with a Terumo Capiox E membrane oxygenator and a Bentley® arterial line filter (25 µm). The roller pump and oxygenator were primed with 2.5 to 3.0l of lactated Ringer's and 1 to 2 units (500 to 1000 ml) of canine whole blood in CPDA-1 solution. The blood was obtained 48 to 72 hours prior to the experiment from a donor animal and refrigerated at 4°C. The animal was systemically heparinized with 300-400 IU·kg⁻¹ of heparin (Organon: porcine intestine®) and additional doses as required, to give an activated clotting time (ACT) ≥ 400 sec (Hemochron 400®). Throughout the experiment the animal had an intravenous infusion of lactated Ringer's at 200-250 ml·hr⁻¹ containing 25 mEq·l⁻¹ of NaHCO₃. This was done to maintain a stable hemoglobin concentration and acid-base state during the experiment (α-stat acid-base management). Norepinephrine (40 µg) was injected into the oxygenator coincident with initiating CPB to minimize the hemodynamic consequences. The animals were randomized to one of two groups: non-pulsatile bypass group; Group NP (n = 6), or computer-controlled bypass group; Group CP.

Group CP (n = 6). Following the initiation of CPB, cooling to 28°C commenced immediately in both groups. Temperature was altered using a Travenol heat exchanger. In both groups of animals, the mean cerebral perfusion pressure (CPP; MAP - mean CSFP) was maintained at greater than 60 mm Hg. Hypothermic non-pulsatile CPB continued

SUBSTITUTE SHEET

2184992

23

for 105 min in Group NP and for 15 min in Group CP while the computer-control was being established (see below) and then for 90 minutes with computer pulsation. At 105 min, rewarming was commenced. In both groups, rewarming to baseline temperature was over 30 min. At 45 min after the start of rewarming, cerebral blood flow (CBF) and blood gas samples were obtained.

In these experiments, for each animal in Group CP, by means of a data acquisition system, a 15 min data file of blood pressure was obtained following induction of anesthesia. Data from a typical modulation file is depicted graphically in Figure 2. A typical output from the computer-controller roller pump relating computer voltage output and the changes in MAP are shown in Figure 3. The data is processed by a computer programmed using Table 1.

Cerebral Blood Flow Measurements: The radioactive microspheres, ultrasonicated in saline, were injected into the arterial cannula, approximately 1 meter proximal to the aortic root, after the PaCO_2 was stable between 35-40 mm Hg. If the PaCO_2 could not be stabilized in this range by adjusting the O_2 flow to the oxygenator, CO_2 was added with a Sechrist® mixer. Approximately 2.5×10^6 microspheres (15 μm diameter) were injected into the arterial cannula. This number of microspheres assured greater than 400 microspheres/sample for accurate blood flow measurement (Heymann et al., 1977). The randomly selected microspheres were labelled with ^{46}Sc , ^{85}Sr , ^{141}Ce , ^{95}Nb , or ^{113}Sn (New England Nuclear). A Harvard pump® withdrew a reference blood sample for determination of organ blood flow (25 ml) from the brachial artery (Compugamma®) after being weighed. The counts/min were converted to regional CBF in $\text{ml} \cdot \text{g}^{-1} \cdot \text{min}^{-1}$ with the use of standard equations.

Total CBF (tCBF) was determined by summing weighted flows to all brain regions and dividing by total brain

SUBSTITUTE SHEET

weight. Similarly, cerebral hemispheric CBF (hCBF) and brain stem CBF (bsCBF) were determined by the summation of weighted flows to the cerebral hemispheres and brain stem, respectively. The CPP was measured as (MAP - mean CSFP) and cerebral metabolic rate for O_2 ($CMRO_2$) as $hCBF \times (Art - SSS O_2 \text{ content})$ in $ml O_2 \cdot g^{-1} \cdot min^{-1}$.

Statistical Analyses: Changes over time for blood gas and hemodynamic variables were evaluated by analysis of variance (ANOVA) for repeated measures. When ANOVA was significant, comparisons were made with the least-squares means test. Data are presented as mean $\pm$ SD.

Example 5

This Example describes the results obtained using the materials and methods described in Example 4.

Temperatures and hemodynamic data are shown in Table 7 below. The temperature did not differ between groups for either the period of hypothermia or following rewarming. In all instances, the nasopharyngeal temperature was able to be increased to $35^\circ C$ within the 30 min time frame without exceeding a temperature gradient of $8^\circ C$ between the heat exchanger and the nasopharyngeal measurement sites. The MAP remained stable over the two temperatures in both groups. A difference in MAP was seen between groups with MAP being greater at both temperatures in Group CP but there was no group x time interaction ($p = 0.0904$). In both groups the CSFP increased with rewarming. The CPP was stable over time, within groups, with no group x time interaction ($p = 0.771$).

The blood gas and blood O_2 content data are shown in Table 8 below. Both groups had similar hemoglobin concentrations during CPB and similar pH and $PaCO_2$. A significant group by time interaction was observed for SSS O_2 content ($p = 0.0005$), SSS O_2 ($p = 0.003$), and art-SSS O_2 content difference ($p = 0.011$). In all instances Group CP remained more stable. In Group NP, there was a

SUBSTITUTE SHEET

2184992

25

significant decrease in SSS O₂ content and SSS O₂ with rewarming, and an increase in the art-SSS O₂ content difference.

The regional CBF and CMRO₂ data are shown in Table 8. In both groups regional CBF increased with rewarming. There was no difference between groups for CBF in any region. Flow:metabolic coupling decreased with rewarming in Group NP. There was no difference in CMRO₂ between groups. In both groups CMRO₂ increased with rewarming.

10 Example 6

This Example discusses the results obtained in Example 5.

Use of computer-controlled CPB roller pump according to the invention, which restores inherent biologic variability, as described in Examples 4 and 5, prevents cerebral deoxygenation during rewarming. The SSS O₂, SSS O₂ content, and the art-SSS O₂ content difference were all stable following rewarming in Group CP. In contrast, in Group NP, the SSS O₂ and SSS O₂ content decreased and the art-SSS O₂ content difference increased with rewarming suggesting cerebral deoxygenation with conventional non-pulsatile CPB (Table 8, Fig. 7). This experimental group had CPB managed similarly to that of patients in a study by Croughwell et al. (Ann. Thorac. Surg. 1994; 58:1702-1708) (α -stat acid-base management, use of arterial line filter and membrane oxygenator, and a similar duration of CPB). In all instances, the cerebral deoxygenation was not as severe as that seen by Croughwell et al. However, the changes seen were similar to those in the clinical scenario, and the more important observation is that computer-controlled CPB prevented these changes with rewarming.

The differences between the two groups does not appear to be a consequence of changes in CBF with computer-controlled CPB. It is of interest that the O₂ content difference was lower during hypothermia in Group

SUBSTITUTE SHEET

NP than in Group CP and became greater with rewarming. This may suggest that there was a difference in distribution of CBF during CPB between the two groups given no difference in tCBF. If so, this finding could
5 imply two flow pathways with CPB (a shunt and a parenchymal flow pathway). The use of computer-controlled CPB which restored inherent biologic rhythms would appear to provide better parenchymal flow, resulting in a greater oxygen content difference during
10 the hypothermic bypass period and an ability for recruitment of the capillary bed in proportion to the requirements of increased metabolic demand with rewarming. With non-pulsatile bypass, the lower O₂ content difference at similar tCBF suggests a greater
15 shunt flow, during the period of hypothermia. With rewarming, parenchymal blood flow appears inadequate to meet the increased metabolic demands of the tissue which results in increased extraction of oxygen and an increased O₂ content difference, lowered SSS O₂ and O₂
20 saturation. These results suggest that the greater deoxygenation seen with conventional CPB may be a consequence of capillary closure to parenchymal beds due to non-pulsatile bypass.

Non-pulsatile perfusion has been demonstrated to
25 increase tissue water and alter vascular properties in other tissues such as the lung. These effects of the abnormal pulsation are unaltered by hypothermia, anesthesia or use of arterial filters. Hence these presumed neural protective interventions may not be
30 helpful and may account for lack of clinical correlates of improved neurologic outcome with their institution. The improvements seen with computer-controlled pulsatile flow provides strong indirect evidence that the microembolic theory inadequately explains why the brain
35 is damaged during CPB.

SUBSTITUTE SHEET

The better cerebral oxygenation occurring with computer-controlled CPB is not likely due to any change in the microembolic load presented to the brain. Pulsation, per se, should not alter the microembolic load, and the CBF was identical in the two groups, at both temperatures suggesting another mechanism independent of microemboli being the causative reason for the difference between groups. Microemboli are felt to be the leading candidates to explain the neurologic and neuropsychologic damage following CPB. The microemboli theory cannot effectively explain the reason for the increased O₂ extraction seen on rewarming. An explanation of parallel flow pathways during CPB (one shunt pathway and one parenchymal flow pathway) which is a consequence of the monotonous regular non-pulsatile blood flow with resultant cerebral capillary bed closure can explain why neural damage occurs despite strategies to decrease cerebral blood flow and thereby decrease the embolic load to the brain. Computer-controlled pulsation creates a more physiologic flow state with improved cerebral oxygenation following rewarming.

SUMMARY OF DISCLOSURE

In summary of this disclosure, the present invention provides computer control of the operation of a cardiopulmonary bypass pump, a lung ventilator or other device which provides simulation of in vivo variability of flow of a biologic fluid to an organ. Modifications are possible within the scope of this invention.

SUBSTITUTE SHEET

Table IProgram Listing for Blood Pump Control

{ \$A+, B-, D+, E+, F-, I+, L+, N+, O-, R+, S+, V+ }

{ \$M 16384, 0, 655360 }

{ Program to control the RPM of a perfusion pump by a variation derived from a file generated by Codas of a pressure recording - before heart bypass has occurred.

Input file is format commented in LoadPeakRates. The Peak pressure and the heart rate period are the only columns used. The minimum of the Peak pressure is removed from each observation to generate a list of Peak Pressure variations from minimum to maximum.

The Peak Pressure variations from 0 to max are scaled from 0 to 5 volts and sent out the D/A channel 0 output. External scaling is used to establish the RPM variation associated with 5 volts output.

The period each pressure variation is held at the output is related to the heart rate period.

The data is scanned forwards and then backwards continuously until a key is struck.

}

Program Purfus;

Uses Crt, Dos;

Const

MAXBUF = 8000; {maximum # of 3 sets of data points to be loaded}

Type

datbuf = array[1..MAXBUF] of single;

datdef = ^datbuf;

statdef = array[1..3] of single;

Var

peak, mean, period : datdef;

peakstat, meanstat, periodstat : statdef; {stats of data loaded}

BaseAddr : word;

{Driver for Das16 D/A channel 0 = BaseAddr + 4, + 5}

SUBSTITUTE SHEET

```
Procedure D2A0Out( x : integer);
var
  LowByte, HighByte : integer;
begin
  LowByte := (xAND $F) * 16;
  HighByte := xDIV 16;
  port[BaseAddr + 4] := LowByte;
  port[BaseAddr + 5] := HighByte; {output not upddated
                                   until High byte sent}
end; { D2A0Out }

Procedure GetHeap;

Procedure Bye;
begin
  writeln('Sorry, not enough memory available for heap!');
  halt(1);
end; { Bye }

begin { GetHeap }
  if MaxAvail < SizeOf(datbuf) then Bye;
  GetMem(peak, SizeOf(datbuf));
  if MaxAvail < SizeOf(datbuf) then Bye;
  GetMem(mean, SizeOf(datbuf));
  if MaxAvail < SizeOf(datbuf) then Bye;
  GetMem(period, SizeOf(datbuf));
end; { GetHeap }

Procedure FreeHeap;
begin
  FreeMem(peak, SizeOf(datbuf));
  FreeMem(mean, SizeOf(datbuf));
  FreeMem(period, SizeOf(datbuf));
end; { FreeHeap }

Function FileExist(fn : string) : boolean;
var
```

SUBSTITUTE SHEET

```
    fp : text;
begin
    Assign(fp, fn);
    {$i} Reset(fp); {$i+}
    if IoResult <> 0 then
        FileExist := false
    else
        begin
            Close(fp);
            FileExist := true;
        end;
end; { FileExist }

Function init : boolean; {test input file & init D/A BaseAddr}
var
    P : PathStr;
    D : DirStr;
    N : NameStr;
    E : ExtStr;
    fp : Text;

begin
    if (ParamCount = 0) or (FileExist(ParamStr(1)) = false) then
        begin
            init := false;
            exit;
        end;

    { Load D/A Configuration file }
    P := ParamStr(0); { full pathname to program running }
    FSplit(P, D, N, E);
    Assign(fp, D + N + '.cfg');
    {$i-} Reset(fp); {$i +}
    if IoResult < > 0 then {file missing, create a default one}

begin
```

SUBSTITUTE SHEET

```

ReWrite(fp);
BaseAddr := 784;
writeln(fp, BaseAddr); (* Al's machine D/A uses & H310 base
address *)
end
else
  readln(fp, BaseAddr);
close(fp);
D2A0Out(0); { send a 0 to D/A to init }
  init := true;
end; { init }
{
      Sample File Format
"Data file: C:\CODAS\MUTCH3"
"Source channel: 6"
"Engineering units: mm Hg"
"Sample rate of channel:" 160
"Cycle type: Peak to Valley"
"Cycles per average": "1
"Valley", " Peak", " Mean", " Sec", " S/N Vly", " S/N Pk"
9.660E+01, 1.335E+02, 1.157E+02, 1.5625E-01, 5.500000E+01,
3.000000E+01
9.635E+01, 1.334E+02, 1.171E+02, 1.4375E-01, 1.160000E+02,
9.300000E+01
}

Procedure LoadPeakRates(fn : string; Var numgot : integer);
var
  fp : text;
  t : string;
  i : integer;
  valley : single;
Function GetReal : single; { Codas file contains[,] so string
bang it)
var
  ierr : integer;
  tmp : string;

```

SUBSTITUTE SHEET

```
x : single;
begin
  while t[1] = ' ' do Delete(t, 1, 1); {leading delete []'s}
  if t[Length(t)] < > ',' then t := t + ','; {force[,] at end of
  str}

  tmp := Copy(t, 1, Pos(',', t) - 1);
  Val(tmp, x, ierr);
  if ierr < > 0 then
  begin
    writeln('an error occurred converting [, tmp, ');
    Close(fp);
    halt(1);
  end;

  Delete(t, 1, Pos(',', t)); {remove value converted from str}
  GetReal := x;
end; { GetReal }

begin { LoadPeakRates }
  Assign(fp, fn);
  {$i-} Reset(fp); {$i+}
  if IoResult < > 0 then halt(1); { file existed previously}

  i := 0; { init line counter }
  numgot := 0;

  while (not Eof(fp)) and (numgot < MAXBUF) do
  begin
    inc(i);
    readln(fp, t);

    if i < 10 then      { list header while loading file }
    begin
      while Pos('"', t) > 0 do { delete []'s for easy reading }
        Delete(t, Pos('"', t), 1);
```

SUBSTITUTE SHEET

2184992

33

```
writeln(t);
end
else
begin
  inc(numgot);
  if numgot <= 10 then writeln(t); { list 1st 10 lines }
  valley := GetReal;
  peak[numgot] := GetReal;
  mean[numgot] := GetReal;
  period[numgot] := GetReal;
  { writeln(valley:13:-4, peak[numgot]:13:-4,
    mean[numgot]:13:-4, period[numgot]:13:-4); }

  end;
end; { while not eof }

Close(fp);

if numgot = MAXBUF then
  writeln('***data buffer only contains first ', MAXBUF,
    ' file entries ***');
end; { LoadPeakRates}

Function Min(x, y : single) : single;
begin
  if x < y then
    Min := x
  else
    Min := y;
end; { Min }

Function Max(x, y : single) : single;
begin
  if x > y then
    Max := x
  else
```

SUBSTITUTE SHEET

```
    Max : = y;  
end;  ( Max )
```

```
Procedure CalStat(Var col : datdef; Var stat : statdef;n :  
integer);  
{ Generaric calculate min, max, and avg of one column }  
var  
    i : integer;  
begin  
    stat[1] : = col^[1];  { init min }  
    stat[2] : = col^[1];  { init max }  
    stat[3] : = 0;        { init summation }  
  
    for i : = 1 to n do  
        begin  
            stat[1] : = Min(stat[1], col^[i]);  
            stat[2] : = Max(stat[2], col^[i]);  
            stat[3] : = stat[3] + col^[i];  
        end;  
    stat[3] : = stat[3] / Int(n);  
end;  ( CalStat )
```

```
Procedure GetStats(n : integer); { cal min, max, & avg of data  
loaded }  
var  
    i : integer;  
begin  
    if n < 2 then exit;  { nothing to do }  
    CalStat(peak, peakstat, n);  
    CalStat(mean, meanstat, n);  
    CalStat(period, periodstat, n);  
  
    { report }  
    writeln;  
    writeln(' ':15, 'Peak':15, 'Mean':15, 'Period:15);  
    writeln('Min':15, peakstat[1]:15:3, meanstat[1]:15:3,
```

SUBSTITUTE SHEET

```
periodstat[1]:15:3, periodstat[1]:15:5);
writeln('Max':15, peakstat[2]:15:3, meanstat[2]:15:3,
periodstat[2]:15:5);
writeln('Avg':15, peakstat [3]:15:3, meanstat[3]:15:3,
periodstat[3]:15:5);
writeln;
end; ( GetStats )

Procedure RemoveMinFromPeak(n : integer);
var
  i : integer;
begin
  for i : = 1 to n do
    peak^[i] := peak^[i] - peakstat[1]; ( remove min from peak
    column )
  peakstat[2] := peakstat[2] - peakstat[1];
  peakstat[1] := 0;
end: ( RemoveMinFromPeak )

Procedure Pause(seconds : single);
Type
  timerec = record
    hr, min, sec, hund : word;
  end;
var
  tnow, tend : timerec;
  t1, t2 : longint;

Function Time2Long( time : timerec) : longint;
var
  t : longint;
  hr, min, sec, hund : longint;
begin
  hr := time.hr;
  min := time.min;
  sec := time.sec;
```

SUBSTITUTE SHEET

2184992

36

```
hund := time.hund;
t := hund + sec * 100 + min * 100 * 60 + hr * 100 * 60 * 60;
Time2Long := t;
end; ( Time2Long )

begin ( Pause )
  with tend do GetTime(hr, min, sec, hund); ( Get Current Time )

  ( add seconds passed in to tend to know stop time )
  with tend do
begin
  hund := hund + Round(Frac(seconds) * 100.0);
  sec := sec + Trunc(Int(seconds));

  if hund > 99 then
begin
  hund := hund - 100;
  sec := sec + 1;
end;

while sec > 59 do
begin
  sec := 60;
  min := min + 1;
end;

while min > 59 do
begin
  min := - 60;
  hr := hr + 1;
end;

while hr > 23 do
  hr := hr - 24;

end; ( with tend )
```

SUBSTITUTE SHEET

2184992

37

```
t2 : = Time2Long(tend);
```

```
repeat
```

```
  with tnow do GetTime(hr, min, sec, hund);  {Get Time now } tnow
```

```
  t1 : = Time2Long(tnow);
```

```
until t1 > = t2;
```

```
end;  { Pause }
```

```
Procedure ListWithPause(n : integer);
```

```
var
```

```
  i, d2a : integer;
```

```
begin
```

```
  for i : = 1 to n do
```

```
begin
```

```
  d2a : = Round(peak^[i] / peakstat[2] * 4095.0);
```

```
  writeln(i:5, peak^[i]:10:1, period^[i]:10:4, d2a:10);
```

```
  Pause(period^[i]);
```

```
end;  { for i : = }
```

```
end;  { ListWithPause }
```

```
Procedure SendToD2A(n : integer);
```

```
var
```

```
  i, d2a : integer;
```

```
  up, done, first : boolean;
```

```
  ch : char;
```

```
begin
```

```
  done : = false;
```

```
  first : = true;  { flag for 1st pass }
```

```
  up : = true  { scan data upwards first }
```

```
  i : =  { data pointer }
```

```
while Not Done do
```

```
begin
```

```
  if (up) and (i = n) then up : = false;  { Up or Down }
```

```
  if (Not up) and (i = 1) then up : = true;
```

SUBSTITUTE SHEET

```
if first then ( 1st pass only )
  first : = false

else ( normal data pointer incr/dec )
begin
  if up then
    Inc(i)
  else
    Dec(i);
end;

d2a : = Round(peak^[i] / peakstat[2] * 4095.0);
D2A0Out(d2a);
Pause(period^[i]);

if Key Pressed then
begin
  ch : = ReadKey;
  if ch = Chr(27) then done : = true;
end;
end; ( while not Done )
end; ( SendToD2A )

var
  n : integer;
  ch : char;

begin ( Main )
  if not init then
  begin
    writeln('Purfus filename.exc');
    halt(1);
  end;

  ClrScr;
  GetHeap;
```

SUBSTITUTE SHEET

```
TextBackGround(1);
TextColor(7);

writeln('loading ', ParamStr(1));
LoadPeakRates(ParamStr(1), n);
writeln('Loaded ', n, ' data points');
GetStats(n);
RemoveMinFromPeak(n); { remove min value from Peak column }
GetStats(n);

writeln('Press any key to continue');
ch := ReadKey;

ClrScr;
writeln('**** Calibration ****');
writeln;
writeln('Set the Pump Rate to desired Blood Pressure');
writeln('Press any key to continue');
ch := ReadKey;
writeln;

writeln('Adjust the Variability Level Control to 0 (CCW) Now!');
writeln('Press any key to continue');
ch := ReadKey;
writeln;

writeln('Set the Variability Level Control to desired MAX Peak
Blood Pressure');
D2A0Out(4095); { set D/A output to 5 volts }
writeln('Press any key to continue');
ch := ReadKey;
D2A0Out(0); { set D/A output to 0 volts }
writeln;

writeln('About to start Pump Control Loop');
writeln('Press any key to continue');
```

SUBSTITUTE SHEET

2184992

40

```
ch : = ReadKey;
writeln;
writeln('Press ESC to terminate Control of Pump Rate');

{ ListWithPause(n); }
  SendToD2A(n);

D2A0Out(0); { reset D/A to 0 on exit }
FreeHeap;
TextBackGround(0);
TextColor(7);
end.
```

SUBSTITUTE SHEET

Vent.Bas - Main of Ventilator Control

21 Nov 93

14:48:59

Microsoft (R) QuickBASIC Compiler Version 4.50

'\$Title:'Vent.Bas - Main of Ventilator Control' \$LineSize:112'

'Vent.Bas

'Chris McLennan

'Oct 11/93

'MenuSys.Bas

DECLARE SUB GetConfig ()

DECLARE SUB Init ()

DECLARE SUB MainMenu ()

DECLARE SUB SetConfig (Mode AS INTEGER)

DECLARE SUB TestGrModes ()

'DasLib.Bas

DECLARE SUB InitDas16 ()

'VentLib.Bas

DECLARE SUB InitVentLib ()

'\$INCLUDE: 'DasLib.Bil'

'----- DasLib.Bil -----

'Das16 Common Area

'Das16 Parameters MUST be in common

CONST A2dBuffSize = 2000

DIM Dio(4) AS INTEGER

DIM A2dBuff(2000) AS INTEGER

DIM A2dChNum(2000) AS INTEGER

COMMON Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

COMMON BiPolar AS INTEGER, DasChan AS INTEGER, VoltFS AS SINGLE

COMMON A2dBuff() AS INTEGER 'Buffer for DasMode4

COMMON A2dChNum() AS INTEGER 'Buffer for A/D channel

info

'for Mode 5 & 6, Quick Basic manages the Data Segment

'\$DYNAMIC

'Dynamic Common is dimensioned AFTER the COMMON specification

CONST A2dBuf2Size = 2000

SUBSTITUTE SHEET

Vent.Bas - Main of Ventilator Control

```

COMMON A2dBuf2() AS INTEGER          'Buffer for A/D data
'$STATIC
'----- End DasLib.Bil -----
'$INCLUDE: 'MenuSys.Bil'
'----- MenuSys.Bil-----
' Constants for best available screen mode
CONST VGA = 12
CONST MCGA = 13
CONST EGA256 = 9
CONST EGA64 = 8
CONST MONO = 10
CONST HERC = 3
CONST CGA = 1
' User-defined type to hold information about the mode
TYPE Config
    Scrn      AS INTEGER
    Colors    AS INTEGER
    Atribs    AS INTEGER
    XPix      AS INTEGER
    YPix      AS INTEGER
    TCOL      AS INTEGER
    TROW      AS INTEGER
END TYPE
'*****Graphics Modes
COMMON VC AS Config, InitRows AS INTEGER, BestMode AS INTEGER,
Available AS STRING
'*****Menu system
'$DYNAMIC
COMMON PopBuf() AS INTEGER, CurPath$, ForGnd AS INTEGER, Bakgnd
AS INTEGER
'$STATIC
'----- End MenuSys.Bil-----
'$INCLUDE: 'VentLib.Bil'
'----- VentLib.Bil -----
'VentLib Dynamic Common declaration

```

SUBSTITUTE SHEET

Vent.Bas - Main of Ventilator Control

```

TYPE VentDat                                'Structure to hold modulation data
    Time AS SINGLE
    Mode AS STRING * 1
    Modulation AS SINGLE
END TYPE
CONST VentArrSize = 2000
COMMON VentArrNum AS INTEGER                '# of items loaded in
VentArr()
COMMON VentDatFn$                            'Data File loaded
'$DYNAMIC
'Dynamic Common is dimensioned AFTER the COMMON specification
COMMON VentArr() AS VentDat                'Buffer for Ventilator
modulation data
'$STATIC
'----- End VentLib.Bi1 -----
'$INCLUDE: 'DasLib.Bi2'
'----- DasLib.Bi2 -----
    DIM A2dBuf2(A2dBuf2Size) AS INTEGER
'----- End DasLib.Bi2 -----
'$INCLUDE: 'MenuSys.Bi2'
'----- MenuSys.Bi2-----
    DIM PopBuf(2002, 1) AS INTEGER 'PopUp Buffer
'----- End MenuSys.Bi2-----
'$INCLUDE: 'VentLib.Bi2'
'----- VentLib.Bi2 -----
    DIM VentArr(VentArrSize) AS VentDat 'Buffer for Ventilator
modulation data
'----- End VentLib.Bi2 -----
    InitDas16                                'General Reset of Das16
    InitVentLib                              'Clear Commons
    'Determine Graphic Card
    GetConfig
    VC.Scrn = BestMode
    SetConfig VC.Scrn                        'loads record VC with graphic mode
parameters

```

SUBSTITUTE SHEET

Vent.Bas - Main of Ventilator Control

' TestGrModes

Init

MainMenu

COLOR 7, 0

END

43933 Bytes Available

41609 Bytes Free

0 Warning Error(s)

0 Severe Error(s)

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

```

'$Title:'DasLib.Bas - Das16 Library' $LineSize:112'
'DasLib
'Das16 Library Interface
'Oct 2/93
'Chris McLennan
'MenuSys.Bas
DECLARE FUNCTION Exist! (FileName$)
DECLARE SUB DAS16 (MODE%, BYVAL dummy%, Flag%)
DECLARE SUB DasMode0 ()
DECLARE SUB DasMode1 (ChLow AS INTEGER, ChHigh AS INTEGER)
DECLARE SUB DasMode2 (NextCh AS INTEGER, ChLow AS INTEGER,
ChHigh AS INTEGER)
DECLARE SUB DasMode3 (A2dData AS INTEGER, A2dCh AS INTEGER)
DECLARE SUB DasMode4 (NumPts AS INTEGER)
DECLARE SUB DasMode5 (NumPts AS INTEGER, Cycle AS INTEGER)
DECLARE SUB DasMode6 (NumPts AS INTEGER, Cycle AS INTEGER)
DECLARE SUB DasMode7 ()
DECLARE SUB DasMode8 (Op AS INTEGER, Status AS INTEGER, Count
AS INTEGER)
DECLARE SUB DasMode9 (NumPts AS INTEGER, StartPt AS INTEGER)
DECLARE SUB DasMode15 (D2aCh AS INTEGER, D2aData AS INTEGER)
DECLARE SUB DasMode16 (D2aDat0 AS INTEGER, D2aDat1 AS INTEGER)
DECLARE SUB DasMode17 (Rate!)
DECLARE SUB DasError (code AS INTEGER)
DECLARE FUNCTION A2dToVolt! (A2dVal AS INTEGER)
DECLARE FUNCTION IntToReal! (x AS INTEGER)
DECLARE FUNCTION RealToInt% (x!)
DECLARE FUNCTION VoltToD2a% (volt!)
'$INCLUDE: 'DasLib.Bil'
'----- DasLib.Bil -----
'Das16 Common Area
'Das16 Parameters MUST be in common
CONST A2dBuffSize = 2000
DIM Dio(4) AS INTEGER
DIM A2dBuff(2000) AS INTEGER

```

SUBSTITUTE SHEET

2184992

DasLib.Bas - Das16 Library

```

DIM A2dChNum(2000) AS INTEGER

COMMON Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER
COMMON BiPolar AS INTEGER, DasChan AS INTEGER, VoltFS AS SINGLE
COMMON A2dBuff() AS INTEGER          'Buffer for DasMode4
COMMON A2dChNum() AS INTEGER          'Buffer for A/D channel
info
'for Mode 5 & 6, Quick Basic manages the Data Segment
'$DYNAMIC
'Dynamic Common is dimensioned AFTER the COMMON specification
CONST A2dBuf2Size = 2000
COMMON A2dBuf2() AS INTEGER          'Buffer for A/D data
'$STATIC
'----- End DasLib.Bi1 -----
'$INCLUDE: 'DasLib.Bi2'
'----- DasLib.Bi2 -----
    DIM A2dBuf2(A2dBuf2Size) AS INTEGER
'----- End DasLib.Bi2 -----
'$Page $SubTitle:'FUNCTION A2dToVolt'
```

SUBSTITUTE SHEET

2184992

47

DasLib.Bas - Das16 Library

FUNCTION A2dToVolt

FUNCTION A2dToVolt (A2dVal AS INTEGER)

'Returns A2d value converted into voltage

SHARED BiPolar AS INTEGER, VoltFS AS SINGLE

IF BiPolar THEN

A2dToVolt = A2dVal / 2048! * VoltFS

ELSE

A2dToVolt = A2dVal / 4096! * VoltFS

END IF

END FUNCTION

'\$Page \$SubTitle:'SUB DasError'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasError

SUB DasError (code AS INTEGER)

SELECT CASE code

CASE 0

PRINT "OK"

CASE 1

PRINT "Driver not initialized"

CASE 2

PRINT "Mode number out of range"

CASE 3

PRINT "Base Address out of range"

CASE 4

PRINT "Interrupt Level out of range"

CASE 5

PRINT "DMA Level out of range"

CASE 6

PRINT "Differential Mux scan limits out of range"

CASE 7

PRINT "Single Ended Mux scan limits out of range"

CASE 8

PRINT "Error Code 8?"

CASE 9

PRINT "A/D Timeout Error - hardware - no EOC"

CASE 10

PRINT "Counter division ratio 0 or 1 in mode 17"

CASE 11

PRINT "Number of conversions <= 0 in modes 4,5, or 6"

CASE 12

PRINT "Counter configuration # out of range in mode 10"

CASE 13

PRINT "Digital output data out of range in mode 13"

CASE 14

PRINT "D/A data out of range in modes 15 or 16"

SUBSTITUTE SHEET

2184992

49

DasLib.Bas - Das16 Library

SUB DasError

```
CASE 15
  PRINT "D/A channel # out of range"
CASE 16
  PRINT "Counter read operation not 0 or 1 in mode 12"
CASE 17
  PRINT "Start conversion negative # in mode 9"
CASE 18
  PRINT "Word count 0 or negative in mode 9"
CASE 19
  PRINT "Trigger mode not 0 or 1 in modes 4,5, or 6"
CASE 20
  PRINT "DMA / Interrupt operation already active in modes
5 or 6"
CASE 21
  PRINT "DMA page wrap around in mode 6"
CASE 22
  PRINT "Hardware failure or installation error [Base
Address?]"
CASE 23
  PRINT "Trigger channel inconsistent with configuration
[mode 19]"
CASE 24
  PRINT "Trigger data out of range [mode 19]"
CASE 25

  PRINT "Slope data not 0 or 1 [mode 19]"
END SELECT

END SUB
```

```
'$Page $SubTitle:'SUB DasMode0'
```

SUBSTITUTE SHEET

2184992

50

DasLib.Bas - Das16 Library

SUB DasMode0

SUB DasMode0

'INIT DAS16 DRIVER

'Use NEW format file Das16Adr.Cfg for Card initialization

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

SHARED BiPolar AS INTEGER, DasChan AS INTEGER, VoltFS AS SINGLE

DIM fp AS INTEGER, Status AS INTEGER

IF NOT Exist("Das16Adr.Cfg") THEN

CLS

PRINT "Das16 Configuration file [Das16Adr.Cfg] not found"

END

END IF

OPEN "Das16Adr.Cfg" FOR INPUT AS #1

INPUT #1, Dio(0) 'get base I/O address

INPUT #1, Dio(1) 'interrupt level

INPUT #1, Dio(2) 'D.M.A. level

INPUT #1, VoltFS 'Voltage Range

CLOSE #1

Flag = 0 'error variable

Md = 0 'mode 0 - initialize

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

END

END IF

'---- Find DASH-16 operating configurations ----

'read status register at base address + 8

Status = INP(Dio(0) + 8)

IF (Status AND &H40) = &H40 THEN ' BIPOLAR OR UNIPOLAR

BiPolar = 0

ELSE

BiPolar = -1

END IF

IF (Status AND &H20) = &H20 THEN

DasChan = 16

SUBSTITUTE SHEET

2184992

51

DasLib.Bas - Das16 Library

SUB DasMode0

ELSE

DasChan = 8

' # CHANNELS

END IF

END SUB

'\$Page \$SubTitle:'SUB DasModel'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasModel

SUB DasModel (ChLow AS INTEGER, ChHigh AS INTEGER)

'SET MULTIPLEXER SCAN LIMITS

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

IF ChLow > ChHigh THEN SWAP ChLow, ChHigh

Dio(0) = ChLow 'LOW LIMIT

Dio(1) = ChHigh 'HIGH LIMIT

Md = 1

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

END

END IF

END SUB

.\$Page \$SubTitle:'SUB DasModel5'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasModel5

SUB DasModel5 (D2aCh AS INTEGER, D2aData AS INTEGER)

'Output Data to one D/A Channel

'If the -5vRev is connected to the D/A Ref In on the DB 37,
'an output from 0 to 5V is generated.

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

Dio(0) = D2aCh 'Channel #[0..1] to send to

Dio(1) = D2aData 'Data to send [0..4095]

Md = 15

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

END

END IF

END SUB

'\$Page \$SubTitle:'SUB DasModel6'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasModel6

SUB DasModel6 (D2aDat0 AS INTEGER, D2aDat1 AS INTEGER)

'Output Data to both D/A Channels

'If the -5vRev is connected to the D/A Ref In on the DB 37,
'an output from 0 to 5V is generated.

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

Dio(0) = D2aDat0 'Channel 0

Dio(1) = D2aDat1 'Channel 1

Md = 16

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

END

END IF

END SUB

'\$Page \$SubTitle:'SUB DasModel7'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasModel7

SUB DasModel7 (Rate)

'COMPUTE SCAN RATE PER CHANNEL FOR A/D

,

' HZ = 1,000,000 / (N1 * N2)

,

' WHERE 2 < N1 & N2 > 65535

,

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

DIM NextCh AS INTEGER, ChLow AS INTEGER, ChHigh AS INTEGER

DIM Sr, HzMin, HzMax, Prod, Seed1, Seed2, NumCh AS INTEGER

'DETERMINE # CHANNELS IN USE

NextCh = 0

ChLow = 0

ChHigh = 0

CALL DasMode2(NextCh, ChLow, ChHigh)

NumCh = ChHigh - ChLow + 1 ' NUMCH

'CAL SCAN RATE PER CHANNEL

Sr = Rate * NumCh

'FIX RATE OVERFLOWS

HzMax = 1000000! / 4

HzMin = 1000000! / (65535! * 65535!)

IF Sr > HzMax THEN Sr = HzMax

IF Sr < HzMin THEN Sr = HzMin

'TRICKY ROUTINE

,

Prod = 1000000! / Sr 'PROD = N1 * N2

Seed1 = 2 'STARTING VALUE OF SEED1

DO

IF Prod / Seed1 <= 65535 THEN EXIT DO 'SEED2 WILL BE
VALID

Seed1 = INT(Seed1 * 2) 'DOUBLE NEXT
GUESS FOR SEED1

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library
SUB DasModel7

```
      IF Seed1 > 65535 THEN Seed1 = 65535      '& TEST FOR  
OVERFLOW  
      LOOP  
      Seed2 = INT(Prod / Seed1)  
  
      'CALL DASH16  
      Md = 17  
  
      Dio(0) = RealToInt(Seed1)  
      Dio(1) = RealToInt(Seed2)  
  
      'PRINT Dio(0), Dio(1), 1000000 / (Seed1 * Seed2)  
  
      CALL DAS16(Md, VARPTR(Dio(0)), Flag)  
      IF Flag <> 0 THEN  
          DasError (Flag)  
  
      END  
      END IF  
END SUB  
  
'$Page $SubTitle:'SUB DasMode2'
```

SUBSTITUTE SHEET

2184992

57

DasLib.Bas - Das16 Library

SUB DasMode2

SUB DasMode2 (NextCh AS INTEGER, ChLow AS INTEGER, ChHigh AS
INTEGER)

'READ CURRENT MULTIPLEXER ADDRESS AND SCAN LIMITS

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

Md = 2

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

END

END IF

NextCh = Dio(0) 'NEXT CHANNEL TO CONVERT

ChLow = Dio(1) 'LOWER SCAN LIMIT

ChHigh = Dio(2) 'UPPER SCAN LIMIT

END SUB

'\$Page \$SubTitle:'SUB DasMode3'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasMode3

SUB DasMode3 (A2dData AS INTEGER, A2dCh AS INTEGER)

'Do one A/D conversion and increment Mux

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

Md = 3

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

END

END IF

A2dData = Dio(0) 'A/D data

A2dCh = Dio(1) 'A/D Channel

END SUB

'\$Page \$SubTitle:'SUB DasMode4'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasMode4

SUB DasMode4 (NumPts AS INTEGER)

'Do NumPts A/D conversions directly to array A2dBuff()

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

SHARED A2dBuff() AS INTEGER

IF NumPts > A2dBuffSize THEN NumPts = A2dBuffSize

Dio(0) = NumPts '# pts to convert

Dio(1) = VARPTR(A2dBuff(0)) 'output array location

Dio(2) = 1 'trigger source, 1=timer,

0=external on IP0

'Note: If the timer is used as a trigger source then holding
input IP0

' low will delay starting conversions.

Md = 4

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

 DasError (Flag)

 END

END IF

END SUB

'\$Page \$SubTitle:'SUB DasMode5'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasMode5

SUB DasMode5 (NumPts AS INTEGER, Cycle AS INTEGER)

'Do NumPts A/D conversions and transfer to memory on INTERRUPT

'This uses Dynamically allocated A2dBuf2 in common data segment

'This runs in the background

'Data is accessed by DasMode9

'Cycle = 0 for One Scan and finish

' = 1 for Continuous Scanning --> DasMode7 to stop

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

SHARED A2dBuf2() AS INTEGER

IF NumPts > A2dBuf2Size THEN NumPts = A2dBuf2Size

Dio(0) = NumPts

'# pts to convert

Dio(1) = VARSEG(A2dBuf2(0))

'output array location

Dio(2) = 1

'trigger source, 1=timer,

0=external on IP0

'Note: If the timer is used as a trigger source then holding
input IP0

' low will delay starting conversions.

Dio(3) = Cycle

'0 = One shot and finish

'1 = Continuous Scanning - Mode7

to stop

Md = 5

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

DasMode7

END

END IF

END SUB

'\$Page \$SubTitle:'SUB DasMode6'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasMode6

SUB DasMode6 (NumPts AS INTEGER, Cycle AS INTEGER)

'Do NumPts A/D conversions and transfer to memory via DMA

'This uses Dynamically allocated A2dBuf2 in common data segment

'This runs in the background

'Data is accessed by DasMode9

'Cycle = 0 for One Scan and finish

' = 1 for Continuous Scanning --> DasMode7 to stop

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

SHARED A2dBuf2() AS INTEGER

IF NumPts > A2dBuf2Size THEN NumPts = A2dBuf2Size

Dio(0) = NumPts '# pts to convert

Dio(1) = VARSEG(A2dBuf2(0)) 'output array location

Dio(2) = 1 'trigger source, 1=timer,

0=external on IP0

'Note: If the timer is used as a trigger source then holding
input IP0

' low will delay starting conversions.

Dio(3) = Cycle '0 = One shot and finish

'1 = Continuous Scanning - Mode7

to stop

Md = 6

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

DasMode7

END

END IF

END SUB

'\$Page \$SubTitle:'SUB DasMode7'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasMode7

SUB DasMode7

'DISABLE DMA/INTERRUPT OPERATION OF A/D

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

Md = 7

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

END

END IF

END SUB

'\$Page \$SubTitle:'SUB DasMode8'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

SUB DasMode8

SUB DasMode8 (Op AS INTEGER, Status AS INTEGER, Count AS INTEGER)

'allows monitor of mode 5,6,18, or 20 Background operation

'Op = 0 none

' = 1 Mode 6 [DMA]

' = 2 Mode 5 [Interrupt]

' = 3 Mode 18 [Interrupt]

' = 4 Mode 20 [Interrupt]

'Status = 0 if Done or 1 in Active

'Count = Current word count [# conversions so far]

' = 0 to (Num Chan * pts per chan) - 1

' ex = 4 channels @ 200 pts per channel = [0..799]

' = seems to indicate the current conversion - which is
not yet done.

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

Md = 8

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

END

END IF

Op = Dio(0) 'Operation type

Status = Dio(1) 'Status of Operation

Count = Dio(2) 'Current word count

END SUB

'\$Page \$SubTitle:'SUB DasMode9'

SUBSTITUTE SHEET

2184992

DasLib.Bas - Das16 Library
SUB DasMode9

SUB DasMode9 (NumPts AS INTEGER, StartPt AS INTEGER)

'Transfer Data from A2dBuf2() to A2dBuff() [from memory to array]

'Channel information transferred into A2dChNum()

'StartPt is starting conversion # - usually zero

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER

SHARED A2dBuf2() AS INTEGER 'input area

SHARED A2dBuff() AS INTEGER 'output area for data

SHARED A2dChNum() AS INTEGER 'output area for channel

#'s

IF NumPts > A2dBuf2Size THEN NumPts = A2dBuf2Size

Dio(0) = NumPts '# pts to convert

Dio(1) = VARSEG(A2dBuf2(0)) 'input array location

Dio(2) = StartPt 'start transferring at
conversion StartPt

Dio(3) = VARPTR(A2dBuff(StartPt)) 'Destination of Data to
transfer

Dio(4) = VARPTR(A2dChNum(StartPt)) 'Destination of chan # Data
to transfer

Md = 9

CALL DAS16(Md, VARPTR(Dio(0)), Flag)

IF Flag <> 0 THEN

DasError (Flag)

DasMode7

END

END IF

END SUB

'\$Page \$SubTitle:'FUNCTION IntToReal'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

FUNCTION IntToReal

FUNCTION IntToReal (x AS INTEGER)

'CONVERT SIGNED INTEGER TO REAL [-32768 <==> 32767]

IF x >= 0 THEN

IntToReal = x

ELSE

IntToReal = x + 65536

END IF

END FUNCTION

'\$Page \$SubTitle:'FUNCTION RealToInt%'

SUBSTITUTE SHEET

DasLib.Bas - Das16 Library

FUNCTION RealToInt%

FUNCTION RealToInt% (x)

'CONVERT REAL TO SIGNED INTEGER [-32768 <==> 32767]

IF x <= 32767 THEN

RealToInt% = x

ELSE

RealToInt% = x - 65536

END IF

END FUNCTION

'\$Page \$SubTitle:'FUNCTION VoltToD2a%'

SUBSTITUTE SHEET

2184992

DasLib.Bas - Das16 Library

FUNCTION VoltToD2a%

FUNCTION VoltToD2a% (volt)

'Calculates D/A value from a voltage from 0 to 5V

'Assumes Das16 5V Reference in use

'Returns integer between 0 & 4095

DIM x

x = volt / 5! * 4096!

IF x < 0 THEN x = 0

IF x > 4095 THEN x = 4095

VoltToD2a% = CINT(x)

END FUNCTION

43949 Bytes Available

37723 Bytes Free

0 Warning Error(s)

0 Severe Error(s)

SUBSTITUTE SHEET

2184992

MenuSys.Bas - QBasic PopDown Menu System

```
'$Title:'MenuSys.Bas      -      QBasic      PopDown      Menu      System'
$LineSize:112'
```

```
' MenuSys.Bas
'
```

```
' Skeleton for PopDown Menu System for QBasic
```

```
' Video Mode is detected
```

```
' Chris McLennan Mar 28/93
```

```
'*****User Interface
```

```
DECLARE FUNCTION CenterStr$ (s$, length AS INTEGER)
```

```
DECLARE FUNCTION Exist! (FileName$)
```

```
DECLARE FUNCTION FUse$ (x AS SINGLE, wide AS INTEGER, dec AS
INTEGER)
```

```
DECLARE FUNCTION LPad$ (s$, length AS INTEGER)
```

```
DECLARE FUNCTION Max (x AS SINGLE, y AS SINGLE)
```

```
DECLARE FUNCTION Min (x AS SINGLE, y AS SINGLE)
```

```
DECLARE FUNCTION RPad$ (s$, length AS INTEGER)
```

```
DECLARE FUNCTION StrTok$ (Srce$, Delim$)
```

```
DECLARE SUB AboutBox ()
```

```
DECLARE SUB BoxIt (x1 AS INTEGER, y1 AS INTEGER, x2 AS INTEGER,
y2 AS INTEGER)
```

```
DECLARE SUB BrList (txt$(), NumLin AS INTEGER)
```

```
DECLARE SUB Init ()
```

```
DECLARE SUB MainMenu ()
```

```
DECLARE SUB MenuBar (Item$(), NumLin AS INTEGER, Choice AS
INTEGER, Ky$, x AS INTEGER, Redraw
AS INTEGER)
```

```
DECLARE SUB MenuExec (cmd AS INTEGER, Done AS INTEGER, RtnCode
AS INTEGER)
```

```
DECLARE SUB MsgBox (txt$(), NumLin AS INTEGER, Ky$, RestoreScrn
AS INTEGER, Pause AS INTEGER)
```

```
DECLARE SUB Pop (x1 AS INTEGER, y1 AS INTEGER, x2 AS INTEGER,
y2 AS INTEGER)
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

```
DECLARE SUB PopUpList (Item$(), NumLin AS INTEGER, Choice AS  
INTEGER, Ky$, SaveScrn AS INTEGE  
R, Center AS INTEGER)
```

```
DECLARE SUB SelectFile (FileName$)
```

```
DECLARE SUB SelectNewPath ()
```

```
DECLARE SUB UnPop ()
```

```
'Graphic Mode Detection
```

```
DECLARE SUB GetConfig ()
```

```
DECLARE SUB SetConfig (mode AS INTEGER)
```

```
DECLARE SUB TestGrModes ()
```

```
'VentLib.Bas
```

```
DECLARE SUB BrPrnData ()
```

```
DECLARE SUB Demo ()
```

```
DECLARE SUB LoadVentDatPrn ()
```

```
DECLARE SUB InitDas16 ()
```

```
DECLARE SUB LoopTest ()
```

```
DECLARE SUB LoopTest5 ()
```

```
DECLARE SUB LoopTest6 ()
```

```
DECLARE SUB OhioFn ()
```

```
DECLARE SUB PlotVentDat ()
```

```
DECLARE SUB VentLoop ()
```

```
'$INCLUDE: 'DasLib.Bil'
```

```
'----- DasLib.Bil -----
```

```
'Das16 Common Area
```

```
'Das16 Parameters MUST be in common
```

```
CONST A2dBuffSize = 2000
```

```
DIM Dio(4) AS INTEGER
```

```
DIM A2dBuff(2000) AS INTEGER
```

```
DIM A2dChNum(2000) AS INTEGER
```

SUBSTITUTE SHEET

2184992

MenuSys.Bas - QBasic PopDown Menu System

```
COMMON Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER
COMMON BiPolar AS INTEGER, DasChan AS INTEGER, VoltFS AS SINGLE
COMMON A2dBuff() AS INTEGER          'Buffer for DasMode4
COMMON A2dChNum() AS INTEGER          'Buffer for A/D channel
info
```

```
'for Mode 5 & 6, Quick Basic manages the Data Segment
```

```
'$DYNAMIC
```

```
'Dynamic Common is dimensioned AFTER the COMMON specification
```

```
CONST A2dBuf2Size = 2000
```

```
COMMON A2dBuf2() AS INTEGER          'Buffer for A/D data
```

```
'$STATIC
```

```
'----- End DasLib.Bil -----
```

```
'$INCLUDE: 'MenuSys.Bil'
```

```
'----- MenuSys.Bil-----
```

```
' Constants for best available screen mode
```

```
CONST VGA = 12
```

```
CONST MCGA = 13
```

```
CONST EGA256 = 9
```

```
CONST EGA64 = 8
```

```
CONST MONO = 10
```

```
CONST HERC = 3
```

```
CONST CGA = 1
```

```
' User-defined type to hold information about the mode
```

```
TYPE Config
```

```
    Scrn      AS INTEGER
```

```
    Colors    AS INTEGER
```

```
    Atribs    AS INTEGER
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

```

    XPix      AS INTEGER
    YPix      AS INTEGER
    TCOL      AS INTEGER
    TROW      AS INTEGER
END TYPE

'*****Graphics Modes
COMMON VC AS Config, InitRows AS INTEGER, BestMode AS INTEGER,
Available AS STRING
'****Menu system
'$DYNAMIC
COMMON PopBuf() AS INTEGER, CurPath$, ForGnd AS INTEGER, Bakgnd
AS INTEGER
'$STATIC

'----- End MenuSys.Bi1-----

'$INCLUDE: 'DasLib.Bi21'
'----- DasLib.Bi2 -----
    DIM A2dBuf2(A2dBuf2Size) AS INTEGER
'----- End DasLib.Bi2 -----

'$INCLUDE: 'MenuSys.Bi2'
'----- MenuSys.Bi2-----
    DIM PopBuf(2002, 1) AS INTEGER 'PopUp Buffer
'----- End MenuSys.Bi2-----

CONST PgmVersion = 1.09

' Menu Data

MainMenuData:
    DATA "File"
    DATA "Demos"
    DATA ""

```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

PullDownMenuData:

```
DATA "About",          1003
DATA "Select path..", 1002
DATA "Open..",         1004
DATA "Import PRN..",   1010
DATA "Browse Data",    1011
DATA "Plot Data",      1012
DATA "Dos Shell",      1001
DATA "Quit",           1000
DATA "", -1
DATA "Init",           2000
DATA "Demo",           2001
DATA "OhioFn",         2002
DATA "LoopTest",       2003
DATA "LoopTest5",      2004
DATA "LoopTest6",      2005
DATA "VentLoop",       2006
DATA "", -1
```

' Error trap to make screen independent

VideoErr:

```
SELECT CASE BestMode      ' Fall through until something works
CASE VGA
    BestMode = MCGA
    Available = "12BD"
CASE MCGA
    BestMode = EGA256
    Available = "12789"
CASE EGA256
    BestMode = CGA
    Available = "12"
CASE CGA
    BestMode = MONO
    Available = "A"
```

SUBSTITUTE SHEET

2184992

73

MenuSys.Bas - QBasic PopDown Menu System

```
CASE MONO
  BestMode = HERC
  Available = "3"
CASE ELSE
  PRINT "Sorry. Graphics not available. Need MsHerc.Com?"
  END
END SELECT
RESUME

' Trap to detect 64K EGA
EGAErr:
  BestMode = EGA64
  Available = "12789"
  RESUME NEXT

' Trap to determine initial number of rows so they can be
restored
RowErr:
  IF InitRows = 50 THEN
    InitRows = 43
    RESUME
  ELSE
    InitRows = 25
    RESUME NEXT
  END IF

'$Page $SubTitle:'SUB AboutBox'
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB AboutBox

SUB AboutBox

REDIM a\$(9)

DIM Ky\$

a\$(0) = " ** Status ** "

a\$(1) = ""

a\$(2) = "Version - " + LPad\$(STR\$(PgmVersion), 7)

a\$(3) = ""

a\$(4) = " Memory Available "

a\$(5) = "String - " + LPad\$(STR\$(FRE("x")), 8)

a\$(6) = "Array - " + LPad\$(STR\$(FRE(-1)), 8)

a\$(7) = "Stack - " + LPad\$(STR\$(FRE(-2)), 8)

a\$(8) = ""

a\$(9) = "Press any key to continue"

MsgBox a\$(), 9, Ky\$, -1, -1

ERASE a\$

END SUB

'\$Page \$SubTitle:'SUB BoxIt'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB BoxIt

SUB BoxIt (x1 AS INTEGER, y1 AS INTEGER, x2 AS INTEGER, y2 AS INTEGER)

DIM Lt\$, Hl\$, Rt\$, Vl\$, Lb\$, Rb\$, top\$, mdl\$, btm\$

DIM r AS INTEGER, wide AS INTEGER

'Line chars

Lt\$ = "┌": Hl\$ = "=": Rt\$ = "┐"

Vl\$ = "│"

Lb\$ = "└": Rb\$ = "┘"

wide = x2 - x1 + 1 - 2

top\$ = Lt\$ + STRING\$(wide, Hl\$) + Rt\$

mdl\$ = Vl\$ + STRING\$(wide, " ") + Vl\$

btm\$ = Lb\$ + STRING\$(wide, Hl\$) + Rb\$

FOR r = y1 TO y2

LOCATE r, x1

IF r = y1 THEN

PRINT top\$;

ELSEIF r = y2 THEN

PRINT btm\$;

ELSE

PRINT mdl\$;

END IF

NEXT r

END SUB

'\$Page \$SubTitle:'SUB BrList'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB BrList

SUB BrList (txt\$(), NumLin AS INTEGER)

'Browse txt\$(NumLin) - press ESC to exit

'txt\$(0) is header on line 1

SHARED ForGnd AS INTEGER, BakGnd AS INTEGER

DIM i AS INTEGER, st AS INTEGER, en AS INTEGER

DIM a\$

IF NumLin = 0 THEN EXIT SUB 'nothing to do

COLOR ForGnd, BakGnd

CLS

'header on line 1

COLOR BakGnd, ForGnd

LOCATE 1, 1

PRINT txt\$(0)

COLOR ForGnd, BakGnd

VIEW PRINT 2 TO 25

'Display list

st = 1

en = Min(CSNG(NumLin), 24)

DO

'update display

CLS

FOR i = st TO en

LOCATE 2 + (i - st), 1

PRINT txt\$(i);

NEXT i

SUBSTITUTE SHEET

2184992

77

MenuSys.Bas - QBasic PopDown Menu System
SUB BrList

```

DO
  a$ = INKEY$
LOOP UNTIL a$ <> ""
SELECT CASE LEN(a$)
  CASE 1
    'ESC to quit
    IF a$ = CHR$(27) THEN EXIT DO
  CASE 2
    SELECT CASE ASC(MID$(a$, 2, 1))
      CASE 71
        'Home
        st = 1
        en = Min(CSNG(NumLin), 24!)
      CASE 79
        'End
        en = NumLin
        st = Max(CSNG(en - 23), 1!)
      CASE 72
        'Up
        st = Max(CSNG(st - 1), 1!)
        en = Min(CSNG(st + 23), CSNG(NumLin))
      CASE 80
        'Down
        en = Min(CSNG(en + 1), CSNG(NumLin))
        st = Max(CSNG(en - 23), 1!)
      CASE 73
        'PgUp
        st = Max(CSNG(st - 24), 1!)
        en = Min(CSNG(st + 23), CSNG(NumLin))
      CASE 81
        'PgDn
        en = Min(CSNG(en + 24), CSNG(NumLin))
        st = Max(CSNG(en - 23), 1!)
      CASE 75
        'Left
      CASE 77
        'Right
    CASE ELSE
      PRINT ASC(MID$(a$, 2, 1))
    END SELECT
  END SELECT
END SELECT

```

SUBSTITUTE SHEET

2184992

78

MenuSys.Bas - QBasic PopDown Menu System

SUB BrList

LOOP

VIEW PRINT

END SUB

'\$Page \$SubTitle:'FUNCTION CenterStr\$'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

FUNCTION CenterStr\$

FUNCTION CenterStr\$ (s\$, length AS INTEGER)

' Adds spaces to start and end of s\$, until s\$ is length chars

IF (LEN(s\$) = length) THEN

CenterStr\$ = s\$

ELSEIF (LEN(s\$) > length) THEN

CenterStr\$ = MID\$(s\$, 1, length)

ELSE

DO

s\$ = " " + s\$ + " "

LOOP UNTIL (LEN(s\$) >= length)

IF (LEN(s\$) > length) THEN s\$ = MID\$(s\$, 1, length)

CenterStr\$ = s\$

END IF

END FUNCTION

'\$Page \$SubTitle:'FUNCTION Exist'

SUBSTITUTE SHEET

2184992

80

MenuSys.Bas - QBasic PopDown Menu System

FUNCTION Exist

FUNCTION Exist (FileName\$)

DIM fp AS INTEGER, size AS LONG

fp = FREEFILE

OPEN FileName\$ FOR RANDOM AS #fp

size = LOF(fp)

CLOSE #fp

IF size = 0 THEN

'File doesn't exist

KILL FileName\$

Exist = 0

ELSE

Exist = -1

END IF

END FUNCTION

'\$Page \$SubTitle:'FUNCTION FUse\$'

SUBSTITUTE SHEET

2184992

MenuSys.Bas - QBasic PopDown Menu System

FUNCTION FUse\$

FUNCTION FUse\$ (x AS SINGLE, wide AS INTEGER, dec AS INTEGER)

'Returns X as a string of length wide, with dec digits after decimal

DIM a\$

DIM i AS INTEGER, j AS INTEGER

a\$ = LTRIM\$(STR\$(x))

'ensure a decimal point exists

IF INSTR(a\$, ".") = 0 THEN a\$ = a\$ + "." + STRING\$(dec, "0")

'add trailing 0's if necessary

i = LEN(a\$)

j = INSTR(a\$, ".")

IF (i - j) < dec THEN a\$ = a\$ + STRING\$(dec - (i - j), "0")

'remove trailing 0's if necessary

i = LEN(a\$)

j = INSTR(a\$, ".")

IF (i - j) > dec THEN a\$ = MID\$(a\$, 1, j + dec)

'pad the correct length

a\$ = LPad\$(a\$, wide)

FUse\$ = a\$

END FUNCTION

'\$Page \$SubTitle:'SUB GetConfig'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB GetConfig

```

'      ===== GetConfig
=====
'      Get the starting number of lines and the video adapter.
'
=====
=====
'
SUB GetConfig STATIC
    SHARED InitRows AS INTEGER, BestMode AS INTEGER, Available AS
    STRING
    ' Assume 50 line display and fall through error
    ' until we get the actual number
    InitRows = 50
    ON ERROR GOTO RowErr
    LOCATE InitRows, 1
    ' Assume best possible screen mode
    BestMode = VGA
    Available = "12789BCD"

    ON ERROR GOTO VideoErr
    ' Fall through error trap until a mode works
    SCREEN BestMode
    ' If EGA, then check pages to see whether more than 64K
    ON ERROR GOTO EGAErr
    IF BestMode = EGA256 THEN SCREEN 8, , 1

    ON ERROR GOTO 0

    ' Reset text mode
    SCREEN 0, , 0
    WIDTH 80, 25

END SUB
'$Page $SubTitle:'SUB Init'
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB Init

SUB Init

```
    SHARED CurPath$, ForGnd AS INTEGER, BakGnd AS INTEGER
    DIM fp AS INTEGER, Choice AS INTEGER
    REDIM Men$(2)
    'Ask Dos for current Path
    SHELL "CD > Dir.Tmp"
    fp = FREEFILE
    OPEN "Dir.Tmp" FOR INPUT AS #fp
        INPUT #fp, CurPath$
    CLOSE #fp
    KILL "Dir.Tmp"
    IF RIGHT$(CurPath$, 1) <> "\" THEN CurPath$ = CurPath$ + "\"
    'Select colour or monochrome
    Men$(0) = "Select your Monitor"
    Men$(1) = "Colour"
    Men$(2) = "Monochrome"
    CLS
    ForGnd = 7
    BakGnd = 0
    ' LOCATE 10, 31
    ' PRINT "Select your Monitor"
    DO
        LOCATE 12, 34
        PopUpList Men$(), 2, Choice, Ky$, -1, 0          'pop
    LOOP UNTIL Ky$ = CHR$(13)

    CLS
    IF Ky$ = CHR$(13) AND Choice = 1 THEN BakGnd = 1      'colour
    COLOR ForGnd, BakGnd
    ERASE Men$
END SUB

'$Page $SubTitle:'FUNCTION LPad$'
```

SUBSTITUTE SHEET

2184992

MenuSys.Bas - QBasic PopDown Menu System

FUNCTION LPad\$

FUNCTION LPad\$ (s\$, length AS INTEGER)

' Adds spaces to s\$ on left, until s\$ is length chars

IF (LEN(s\$) = length) THEN

LPad\$ = s\$

ELSEIF (LEN(s\$) > length) THEN

LPad\$ = MID\$(s\$, 1, length)

ELSE

DO

s\$ = " " + s\$

LOOP UNTIL (LEN(s\$) = length)

LPad\$ = s\$

END IF

END FUNCTION

'\$Page \$SubTitle:'SUB MainMenu'

SUBSTITUTE SHEET

2184992

MenuSys.Bas - QBasic PopDown Menu System

SUB MainMenu

SUB MainMenu

SHARED CurPath\$, ForGnd AS INTEGER, BakGnd AS INTEGER

DIM i AS INTEGER, Bchoice AS INTEGER, numpickM AS INTEGER

DIM Done AS INTEGER, tmp AS INTEGER, barX AS INTEGER

DIM DataSet AS INTEGER, Redraw AS INTEGER, RtnCode AS INTEGER

DIM numpick AS INTEGER, Choice AS INTEGER

DIM m\$, Mky\$, Ky\$, CmdLst\$

RESTORE MainMenuData 'determine # picks

numpickM = 0

DO

READ m\$

numpickM = numpickM + 1

LOOP UNTIL m\$ = ""

numpickM = numpickM - 1

REDIM MBar\$(numpickM)

RESTORE MainMenuData 'Load the menu data and cmd table

FOR i = 1 TO numpickM

READ MBar\$(i)

NEXT i

'Init pull down valid command table

'enter, esc, left & right arrow

 CmdLst\$ = CHR\$(13) + CHR\$(27) + (CHR\$(0) + CHR\$(75)) + (CHR\$(0)
+ CHR\$(77))

Done = 0

Redraw = 0

DO

DO

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB MainMenu

MenuBar MBar\$(), numpickM, Bchoice, Mky\$, barX, Redraw
LOOP UNTIL Mky\$ = CHR\$(13)

'count # picks in the selected pull down menu
RESTORE PullDownMenuData

DataSet = 0

DO WHILE DataSet < Bchoice

numpick = 0 'scan group of pull down menu
data

DO

READ m\$, tmp

numpick = numpick + 1

LOOP UNTIL tmp = -1

numpick = numpick - 1 'correct # items in group

DataSet = DataSet + 1 'correct ptr to pull down group

LOOP

REDIM Menu\$(numpick), cmd(numpick) AS INTEGER 'allocate
memory

'Load the selected pull down menu

RESTORE PullDownMenuData

DataSet = 0

'scan group of pull down menu data, until at start of desired
one

DO UNTIL DataSet = Bchoice - 1

DO

SUBSTITUTE SHEET

87 2184992

MenuSys.Bas - QBasic PopDown Menu System

SUB MainMenu

```
    READ m$, tmp
    LOOP UNTIL tmp = -1

    DataSet = DataSet + 1      'correct ptr to pull down group

LOOP

FOR i = 1 TO numpick          'load pull down data
    READ Menu$(i), cmd(i)
NEXT i

VIEW PRINT 2 TO 25            'clear screen under menu bar
COLOR ForGnd, 0
CLS 2
COLOR ForGnd, BakGnd
VIEW PRINT

DO
    LOCATE 2, barX
    PopUpList Menu$(), numpick, Choice, Ky$, 0, 0      'no
pop, no center
    LOOP UNTIL INSTR(CmdLst$, Ky$) > 0

SELECT CASE Ky$

CASE CHR$(0) + CHR$(77)      'right arrow
    Bchoice = Bchoice + 1
    IF Bchoice > numpickM THEN Bchoice = 1
    Redraw = -1              'unpop pulldown & move to adjacent pull
down

CASE CHR$(0) + CHR$(75)      'left arrow
    Bchoice = Bchoice - 1
```

SUBSTITUTE SHEET

88 2184992

MenuSys.Bas - QBasic PopDown Menu System

SUB MainMenu

IF Bchoice < 1 THEN Bchoice = numpickM

Redraw = -1 'unpop pulldown & move to adjacent pull
down

CASE CHR\$(27)

'ESC

Redraw = 0

CASE CHR\$(13)

'Enter

MenuExec cmd(Choice), Done, RtnCode

Redraw = -1

END SELECT

LOOP UNTIL Done

ERASE Menu\$, cmd

'deallocate pull down data

ERASE MBar\$

END SUB

'\$Page \$SubTitle:'FUNCTION Max'

SUBSTITUTE SHEET

2184992

MenuSys.Bas - QBasic PopDown Menu System

FUNCTION Max

FUNCTION Max (x AS SINGLE, y AS SINGLE)

IF x > y THEN

Max = x

ELSE

Max = y

END IF

END FUNCTION

'\$Page \$SubTitle:'SUB MenuBar'

SUBSTITUTE SHEET

2184992

MenuSys.Bas - QBasic PopDown Menu System

SUB MenuBar

SUB MenuBar (Item\$(), NumLin AS INTEGER, Choice AS INTEGER, Ky\$,
x AS INTEGER, Redraw AS INTEGE
R)

' Horizontal Menu Bar
' Selected by down arrow or enter
' if Redraw true, the menu bar is only updated
' Item\$(0) is ignored

SHARED ForGnd AS INTEGER, BakGnd AS INTEGER

DIM maxLen AS INTEGER, i AS INTEGER
DIM hot AS INTEGER, oldhot AS INTEGER, Done AS INTEGER
DIM mask\$, cmd\$, CmdTbl\$

maxLen = 0
FOR i = 1 TO NumLin 'get max width
IF LEN(Item\$(i)) > maxLen THEN maxLen = LEN(Item\$(i))
NEXT i
IF maxLen + NumLin * 2 > 80 THEN STOP

'generate table for column positions for each pick

REDIM whereX(NumLin) AS INTEGER

whereX(1) = 1
FOR i = 2 TO NumLin
whereX(i) = whereX(i - 1) + LEN(Item\$(i - 1)) + 2
NEXT i

VIEW PRINT 1 TO 1
COLOR BakGnd, ForGnd
CLS 2

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB MenuBar

```
'init command table for 1st letter in Item$
CmdTbl$ = ""
FOR i = 1 TO NumLin
    CmdTbl$ = CmdTbl$ + UCASE$(LEFT$(Item$(i), 1))
NEXT i
CmdTbl$ = CmdTbl$ + CmdTbl$           'double list to wrap
around

IF Choice >= 1 AND Choice <= NumLin THEN
    hot = Choice

ELSE
    hot = 1
END IF
oldhot = hot

'Draw MenuBar
FOR i = 1 TO NumLin
    LOCATE 1, whereX(i)
    mask$ = "\" + STRING$(LEN(Item$(i)), " ") + "\"
    IF i = hot THEN COLOR ForGnd, BakGnd
    PRINT USING mask$; " " + Item$(i) + " ";
    IF i = hot THEN COLOR BakGnd, ForGnd

NEXT i

'if Redraw, were done
IF Redraw THEN
    Ky$ = CHR$(13)
    x = whereX(hot)
    VIEW PRINT
    ERASE whereX
    EXIT SUB
END IF
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System
SUB MenuBar

Done = 0

DO

IF hot <> oldhot THEN

'update menu bar

'deselect oldhot

LOCATE 1, whereX(oldhot)

mask\$ = "\" + STRING\$(LEN(Item\$(oldhot)), " ") + "\"

PRINT USING mask\$; " " + Item\$(oldhot) + " ";

'select hot

LOCATE 1, whereX(hot)

mask\$ = "\" + STRING\$(LEN(Item\$(hot)), " ") + "\"

COLOR ForGnd, BakGnd

PRINT USING mask\$; " " + Item\$(hot) + " ";

COLOR BakGnd, ForGnd

END IF

Choice = hot

oldhot = hot

DO

cmd\$ = UCASE\$(INKEY\$)

LOOP UNTIL cmd\$ <> ""

SELECT CASE LEN(cmd\$)

CASE 1 'one char

Ky\$ = cmd\$ 'save key pressed, to return it

SELECT CASE ASC(cmd\$)

CASE 27 'Esc

Done = -1

CASE 13 'Enter

Done = -1

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB MenuBar

```
CASE ELSE                '1st letter of command?
    i = INSTR(hot - 1 + 2, CmdTbl$, UCASE$(cmd$))
    IF i > 0 THEN
        hot = (i - 1) + 1
        IF hot > NumLin THEN hot = hot - NumLin
    END IF
END SELECT

CASE 2                    'extended char
    SELECT CASE ASC(MID$(cmd$, 2, 1))
        CASE 71            'home
            hot = 1
        CASE 79            'end
            hot = NumLin
        CASE 80            'dn
            Ky$ = CHR$(13)    'dn = enter
            Done = -1
        CASE 77            'right
            hot = hot + 1
            IF hot > NumLin THEN hot = 1
        CASE 75            'left
            hot = hot - 1
            IF hot < 1 THEN hot = NumLin
    END SELECT
END SELECT
LOOP UNTIL Done
'return left x position of item selected on menu bar
x = whereX(hot)
VIEW PRINT
ERASE whereX
END SUB

'$Page $SubTitle:'SUB MenuExec'
```

SUBSTITUTE SHEET

2184992

94

MenuSys.Bas - QBasic PopDown Menu System

SUB MenuExec

```
SUB MenuExec (cmd AS INTEGER, Done AS INTEGER, RtnCode AS
INTEGER) STATIC
```

```
' cmd          = command to execute
' done          = if set to -1, causes menu to terminate program
' RtnCode       = error return of command (-1 if fail)
```

```
SHARED CurPath$, ForGnd AS INTEGER, BakGnd AS INTEGER
DIM FileName$
```

```
SELECT CASE cmd
```

```
    CASE 1000                                'quit
        Done = -1
```

```
    CASE 1001                                'shell
        COLOR 7, 0
        CLS
        PRINT "Type Exit to return to program"
        SHELL
        CLS
        COLOR ForGnd, BakGnd
```

```
    CASE 1002                                'new path
        SelectNewPath
```

```
    CASE 1003                                'about
        AboutBox
```

```
    CASE 1004                                'Open
        SelectFile FileName$
```

```
    CASE 1010                                'Import PRN..
        LoadVentDatPrn
```

SUBSTITUTE SHEET

2184992

95

MenuSys.Bas - QBasic PopDown Menu System

SUB MenuExec

CASE 1011

'Browse Data

BrPrnData

CASE 1012

'Plot Data

PlotVentDat

CASE 2000

'Init Das16

InitDas16

CASE 2001

'Demo of Acquisition

primitives

Demo

CASE 2002

'Demo of Ohio Control

Regression Analysis

OhioFn

CASE 2003

'Slow LoopTest

LoopTest

CASE 2004

'Medium speed LoopTest

LoopTest5

CASE 2005

'High Speed LoopTest

LoopTest6

CASE 2006

'Interrupt LoopTest

VentLoop

END SELECT

END SUB

'\$Page \$SubTitle:'FUNCTION Min'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

FUNCTION Min

FUNCTION Min (x AS SINGLE, y AS SINGLE)

IF x < y THEN

Min = x

ELSE

Min = y

END IF

END FUNCTION

'\$Page \$SubTitle:'SUB MsgBox'

SUBSTITUTE SHEET

2184992

97

MenuSys.Bas - QBasic PopDown Menu System

SUB MsgBox

```
SUB MsgBox (txt$(), NumLin AS INTEGER, Ky$, RestoreScrn AS
INTEGER, Pause AS INTEGER)
```

```
' if RestoreScrn NOT true, user must call UnPop
' if Pause true, inkey$ returns Ky$
' Title of box in Txt$(0)
```

```
SHARED CurPath$, ForGnd AS INTEGER, BakGnd AS INTEGER
DIM x1 AS INTEGER, y1 AS INTEGER, x2 AS INTEGER, y2 AS INTEGER
DIM i AS INTEGER, maxLen AS INTEGER
```

```
'calcute area of centered box
```

```
maxLen = 0
```

```
FOR i = 0 TO NumLin
```

```
    IF LEN(txt$(i)) > maxLen THEN maxLen = LEN(txt$(i))
```

```
NEXT i
```

```
IF maxLen > 78 THEN maxLen = 78
```

```
x1 = 41 - (maxLen + 2) \ 2
```

```
x2 = x1 + maxLen + 1
```

```
y1 = 12 - (NumLin + 2) \ 2
```

```
y2 = y1 + (NumLin + 1)
```

```
Pop x1, y1, x2, y2
```

```
BoxIt x1, y1, x2, y2
```

```
'Title - Txt$(0) not blank
```

```
IF txt$(0) <> "" THEN
```

```
    LOCATE y1, (x1 + (x2 - x1 + 1) \ 2) - LEN(txt$(0)) \ 2
```

```
    PRINT txt$(0);
```

```
END IF
```

```
FOR i = 1 TO NumLin
```

```
    LOCATE y1 + i, x1 + 1
```

SUBSTITUTE SHEET

98 2184992

MenuSys.Bas - QBasic PopDown Menu System

SUB MsgBox

```
PRINT CenterStr$(txt$(i), maxLen);  
NEXT i
```

```
IF Pause THEN
```

```
DO
```

```
    Ky$ = UCASE$(INKEY$)
```

```
    LOOP UNTIL Ky$ <> ""
```

```
END IF
```

```
IF RestoreScrn THEN UnPop
```

```
END SUB
```

```
'$Page $SubTitle:'SUB Pop'
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB Pop

SUB Pop (x1 AS INTEGER, y1 AS INTEGER, x2 AS INTEGER, y2 AS
INTEGER)

SHARED PopBuf() AS INTEGER

DIM r AS INTEGER, c AS INTEGER, i AS INTEGER

PopBuf(1, 0) = x1

PopBuf(1, 1) = y1

PopBuf(2, 0) = x2

PopBuf(2, 1) = y2

i = 3

FOR r = y1 TO y2

FOR c = x1 TO x2

PopBuf(i, 0) = SCREEN(r, c, 0)

PopBuf(i, 1) = SCREEN(r, c, 1)

i = i + 1

NEXT c

NEXT r

END SUB

'\$Page \$SubTitle:'SUB PopUpList'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB PopUpList

SUB PopUpList (Item\$(), NumLin AS INTEGER, Choice AS INTEGER,
Ky\$, SaveScrn AS INTEGER, Center
AS INTEGER)

' if SaveScrn TRUE, screen is popped/unpopped
' if Center TRUE, LIST is centered
' Left & Right arrow exit, with Ky\$ set to [0 75], [0 77]
respectivly
' if Choice in range, list is scrolled to selected item

SHARED ForGnd AS INTEGER, BakGnd AS INTEGER

DIM whereX AS INTEGER, whereY AS INTEGER, maxLen AS INTEGER,
i AS INTEGER

DIM x1 AS INTEGER, y1 AS INTEGER, x2 AS INTEGER

DIM y2 AS INTEGER, ps AS INTEGER, pe AS INTEGER, hot AS INTEGER

DIM oldhot AS INTEGER, Done AS INTEGER, scroll AS INTEGER

DIM mask\$, CmdTbl\$

maxLen = 0

whereY = CSRLIN 'save x,y position @ calling

whereX = POS(1)

FOR i = 0 TO NumLin 'get max width

IF LEN(Item\$(i)) > maxLen THEN maxLen = LEN(Item\$(i))

NEXT i

IF maxLen > 78 THEN STOP

'choose popup area

IF NOT Center THEN

y1 = whereY

IF whereY + NumLin + 1 > 24 THEN 'determine lower screen
limits

y2 = 24

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB PopUpList

```

ELSE
    y2 = whereY + NumLin + 1
END IF

x1 = whereX                                'cal popup hort region
x2 = whereX + maxLen + 1
IF x2 > 80 THEN                             'force col 80 = right edge
    x2 = 80
    x1 = 80 - maxLen - 1
END IF

ELSE                                         'Center popup

    x1 = 40 - (maxLen + 2) \ 2
    x2 = x1 + maxLen + 1
    y1 = 12 - (NumLin + 2) \ 2
    IF y1 < 2 THEN y1 = 2
    y2 = y1 + (NumLin - 1 + 2)
    IF y2 > 24 THEN y2 = 24
END IF

COLOR ForGnd, BakGnd
IF SaveScrn THEN Pop x1, y1, x2, y2      'save region on screen
BoxIt x1, y1, x2, y2                     'box it

'Title - if not blank
IF Item$(0) <> "" THEN
    LOCATE y1, (x1 + (x2 - x1 + 1) \ 2) - LEN(Item$(0)) \ 2
    PRINT Item$(0);
END IF

mask$ = "\" + STRING$(maxLen - 2, " ") + "\"

ps = 1

```

SUBSTITUTE SHEET

2184992

102

MenuSys.Bas - QBasic PopDown Menu System

SUB PopUpList

```
pe = ps + (y2 - y1 - 2)
```

```
'if value of Choice called with is in range - use it
```

```
IF Choice >= ps AND Choice <= NumLin THEN
```

```
    hot = Choice
```

```
    DO WHILE hot > pe
```

```
        ps = ps + 1
```

```
        pe = pe + 1
```

```
    LOOP
```

```
ELSE
```

```
    hot = ps
```

```
END IF
```

```
oldhot = hot
```

```
Done = 0
```

```
scroll = 0
```

```
'init command table for 1st letter in Item$
```

```
CmdTbl$ = ""
```

```
FOR i = 1 TO NumLin
```

```
    CmdTbl$ = CmdTbl$ + UCASE$(LEFT$(Item$(i), 1))
```

```
NEXT i
```

```
CmdTbl$ = CmdTbl$ + CmdTbl$           'double list to wrap
around
```

```
FOR i = ps TO pe                       'draw list
```

```
    LOCATE y1 + 1 + i - ps, x1 + 1
```

```
    IF i = hot THEN COLOR BakGnd, ForGnd
```

```
    PRINT USING mask$; Item$(i);
```

```
    IF i = hot THEN COLOR ForGnd, BakGnd
```

```
NEXT i
```

```
DO
```

```
    'Display picks
```

```
    IF hot >= ps AND hot <= pe THEN
```

```
        'pick on screen
```

```
        LOCATE y1 + 1 + oldhot - ps, x1 + 1
```

```
        'deselect old
```

```
        PRINT USING mask$; Item$(oldhot);
```

SUBSTITUTE SHEET

103

2184992

MenuSys.Bas - QBasic PopDown Menu System

SUB PopUpList

```

        LOCATE y1 + 1 + hot - ps, x1 + 1           'select new
        COLOR BakGnd, ForGnd
        PRINT USING mask$; Item$(hot);
        COLOR ForGnd, BakGnd

    ELSE                                           'scroll or moved
by 1 line

        IF scroll THEN                             'scroll list

            scroll = 0

            IF hot < ps THEN                       'scroll up
                ps = ps - (y2 - y1 - 1)
                IF ps < 1 THEN ps = 1
                pe = ps + (y2 - y1 - 2)

            ELSE                                   'scroll down
                pe = pe + (y2 - y1 - 1)
                IF pe > NumLin THEN pe = NumLin
                ps = pe - (y2 - y1 - 2)
            END IF

        ELSEIF hot < ps THEN                       'scroll up 1
            ps = hot
            pe = ps + (y2 - y1 - 2)

        ELSE                                       'scroll down 1
            pe = hot
            ps = pe - (y2 - y1 - 2)

            IF ps < 1 THEN
                ps = 1
                pe = ps + (y2 - y1 - 2)
            END IF

```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB PopUpList

END IF

FOR i = ps TO pe 'draw list

LOCATE y1 + 1 + i - ps, x1 + 1

IF i = hot THEN COLOR BakGnd, ForGnd

PRINT USING mask\$; Item\$(i);

IF i = hot THEN COLOR ForGnd, BakGnd

NEXT i

END IF

Choice = hot

oldhot = hot

DO

cmd\$ = UCASE\$(INKEY\$)

LOOP UNTIL cmd\$ <> ""

SELECT CASE LEN(cmd\$)

CASE 1 'one char

Ky\$ = cmd\$ 'save key pressed, to return it

SELECT CASE ASC(cmd\$)

CASE 27 'Esc

Done = -1

CASE 13 'Enter

Done = -1

CASE ELSE '1st letter of command?

i = INSTR(hot - 1 + 2, CmdTbl\$, UCASE\$(cmd\$))

IF i > 0 THEN

hot = (i - 1) + 1

IF hot > NumLin THEN hot = hot - NumLin

END IF

END SELECT

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB PopUpList

```
CASE 2          'extended char
  SELECT CASE ASC(MID$(cmd$, 2, 1))
    CASE 71      'home
      hot = 1
    CASE 79      'end
      hot = NumLin
    CASE 73      'pg up
      hot = hot - (y2 - y1 - 1)
      IF hot < 1 THEN hot = 1
      scroll = -1
    CASE 81      'pg dn
      hot = hot + (y2 - y1 - 1)
      IF hot > NumLin THEN hot = NumLin
      scroll = -1
    CASE 72      'up
      hot = hot - 1
      IF hot < 1 THEN hot = NumLin
    CASE 80      'dn
      hot = hot + 1
      IF hot > NumLin THEN hot = 1
    CASE 77      'right
      Ky$ = cmd$
      Done = -1
    CASE 75      'left
      Ky$ = cmd$
      Done = -1
  END SELECT
END SELECT
```

LOOP UNTIL Done

IF SaveScrn THEN UnPop

END SUB

'\$Page \$SubTitle:'FUNCTION RPad\$'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

FUNCTION RPad\$

FUNCTION RPad\$ (s\$, length AS INTEGER)

' Adds spaces to s\$ on right, until s\$ is length chars

IF (LEN(s\$) = length) THEN

 RPad\$ = s\$

ELSEIF (LEN(s\$) > length) THEN

 RPad\$ = MID\$(s\$, 1, length)

ELSE

 DO

 s\$ = s\$ + " "

 LOOP UNTIL (LEN(s\$) = length)

 RPad\$ = s\$

END IF

END FUNCTION

'\$Page \$SubTitle:'SUB SelectFile'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB SelectFile

SUB SelectFile (FileName\$)

 SHARED CurPath\$, ForGnd AS INTEGER, BakGnd AS INTEGER

 DIM fp AS INTEGER, nfil AS INTEGER, Choice AS INTEGER, size AS
LONG

 DIM FilSpec\$, Ky\$

 FileName\$ = ""

 REDIM Msg\$(3)

 Msg\$(0) = "Select File"

 Msg\$(1) = ""

 Msg\$(2) = "Scanning " + CurPath\$

 Msg\$(3) = ""

 MsgBox Msg\$(), 3, Ky\$, 0, 0 'no unpop, no pause

 ERASE Msg\$

 SHELL "Dir " + CurPath\$ + " *.* /ON /A-D /B > Dir.Tmp"

 fp = FREEFILE

 OPEN "Dir.Tmp" FOR RANDOM AS #fp LEN = 1

 size = LOF(fp)

 CLOSE #fp

 IF size = 0 THEN

 UnPop

 'remove MsgBox message

 REDIM Msg\$(7)

 Msg\$(0) = "Select File"

 Msg\$(1) = ""

 Msg\$(2) = "----- WARNING -----"

 Msg\$(3) = ""

SUBSTITUTE SHEET

108 2184992

MenuSys.Bas - QBasic PopDown Menu System

SUB SelectFile

```
Msg$(4) = CurPath$
Msg$(5) = "contained NO files"
Msg$(6) = ""
Msg$(7) = "Press any key to continue"
MsgBox Msg$(), 7, Ky$, -1, -1    'unpop, pause
ERASE Msg$
```

KILL "Dir.Tmp"

EXIT SUB

END IF

```
nfil = 0
OPEN "Dir.Tmp" FOR INPUT AS #fp
WHILE NOT EOF(fp)
    LINE INPUT #fp, FilSpec$
    FilSpec$ = LTRIM$(RTRIM$(FilSpec$))
    IF LEN(FilSpec$) > 0 AND FilSpec$ <> "DIR.TMP" THEN nfil =
nfil + 1
WEND
CLOSE #fp
```

IF nfil = 0 THEN

UnPop 'remove MsgBox message

```
REDIM Msg$(7)
Msg$(0) = "Select File"
Msg$(1) = ""
Msg$(2) = "----- WARNING -----"
Msg$(3) = ""
Msg$(4) = CurPath$
Msg$(5) = "contained NO files"
Msg$(6) = ""
Msg$(7) = "Press any key to continue"
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB SelectFile

```
MsgBox Msg$(), 7, Ky$, -1, -1      'unpop, pause
ERASE Msg$

KILL "Dir.Tmp"
EXIT SUB
END IF

REDIM Fil$(nfil)

fp = FREEFILE
OPEN "Dir.Tmp" FOR INPUT AS #fp
nfil = 1
WHILE NOT EOF(fp)
    LINE INPUT #fp, Fil$(nfil)
    Fil$(nfil) = LTRIM$(RTRIM$(Fil$(nfil)))
    IF LEN(Fil$(nfil)) > 0 AND Fil$(nfil) <> "DIR.TMP" THEN nfil
= nfil + 1
WEND
nfil = nfil - 1
CLOSE #fp

KILL "Dir.Tmp"
UnPop                                'remove MsgBox message
Fil$(0) = CurPath$ + ".*"
DO
    PopUpList Fil$(), nfil, Choice, Ky$, -1, -1      'pop,
center
    LOOP UNTIL Ky$ = CHR$(13) OR Ky$ = CHR$(27)

    IF Ky$ = CHR$(13) THEN FileName$ = Fil$(Choice)

    ERASE Fil$
END SUB
'$Page $SubTitle:'SUB SelectNewPath'
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB SelectNewPath

SUB SelectNewPath

 SHARED CurPath\$, ForGnd AS INTEGER, BakGnd AS INTEGER

 DIM fp AS INTEGER, ndir AS INTEGER, Choice AS INTEGER

 DIM Drv\$, DirSpec\$, Ky\$

 REDIM Msg\$(4)

 Msg\$(0) = "Select New Path"

 Msg\$(1) = ""

 Msg\$(2) = "Current Path: " + CurPath\$

 Msg\$(3) = ""

 Msg\$(4) = "Drive to Scan [A..Z, ESC] = "

 DO

 MsgBox Msg\$(), 4, Drv\$, -1, -1

 LOOP UNTIL INSTR("ABCDEFGHIJKLMNOPQRSTUVWXYZ" + CHR\$(13) +
CHR\$(27), Drv\$) > 0

 ERASE Msg\$

 IF Drv\$ = CHR\$(13) THEN

 Drv\$ = LEFT\$(CurPath\$, 1)

 ELSEIF Drv\$ = CHR\$(27) THEN

 EXIT SUB

 END IF

 REDIM Msg\$(3)

 Msg\$(0) = "Select New Path"

 Msg\$(1) = ""

 Msg\$(2) = "Scanning Directory on " + Drv\$ + ":"

 Msg\$(3) = ""

 MsgBox Msg\$(), 3, Ky\$, 0, 0 'no unpop, no pause

 ERASE Msg\$

SUBSTITUTE SHEET

111

MenuSys.Bas - QBasic PopDown Menu System

SUB SelectNewPath

SHELL "Dir " + Drv\$ + ":\ /AD /S /B | Sort > Dir.Tmp"

fp = FREEFILE

ndir = 1

OPEN "Dir.Tmp" FOR INPUT AS #fp

WHILE NOT EOF(fp)

LINE INPUT #fp, DirSpec\$

ndir = ndir + 1

WEND

CLOSE #fp

IF ndir = 2 AND LEN(DirSpec\$) = 0 THEN

'Drive had no

sub-dirs

ndir = 1

END IF

REDIM Fdir\$(ndir)

Fdir\$(1) = Drv\$ + ":\"

'Root Dir

IF ndir > 1 THEN

fp = FREEFILE

OPEN "Dir.Tmp" FOR INPUT AS #fp

ndir = 2

WHILE NOT EOF(fp)

LINE INPUT #fp, Fdir\$(ndir)

ndir = ndir + 1

WEND

ndir = ndir - 1

CLOSE #fp

END IF

KILL "Dir.Tmp"

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB SelectNewPath

UnPop

'remove MsgBox message

IF ndir = 1 THEN

REDIM Msg\$(6)

Msg\$(0) = "Select New Path"

Msg\$(1) = ""

Msg\$(2) = "----- WARNING -----"

Msg\$(3) = ""

Msg\$(4) = Drv\$ + ": had NO directories"

Msg\$(5) = ""

Msg\$(6) = "Press any key to continue"

MsgBox Msg\$(), 6, Ky\$, -1, -1 'unpop, pause

ERASE Msg\$

END IF

Fdir\$(0) = "Select New Path"

DO

PopUpList Fdir\$(), ndir, Choice, Ky\$, -1, -1 'pop,
center

LOOP UNTIL Ky\$ = CHR\$(13) OR Ky\$ = CHR\$(27)

IF Ky\$ = CHR\$(13) THEN

CurPath\$ = Fdir\$(Choice)

IF RIGHT\$(CurPath\$, 1) <> "\" THEN CurPath\$ = CurPath\$ + "\"

END IF

ERASE Fdir\$

END SUB

'\$Page \$SubTitle:'SUB SetConfig'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB SetConfig

```
' ===== SetConfig
=====
' Sets the correct values for each field of the VC variable.
They
' vary depending on Mode and on the current configuration.
'
=====
=====
```

SUB SetConfig (mode AS INTEGER) STATIC

SHARED VC AS Config, BestMode AS INTEGER

SELECT CASE mode

CASE 1 ' Four-color graphics for CGA, EGA, VGA, and MCGA

IF BestMode = CGA OR BestMode = MCGA THEN

VC.Colors = 0

ELSE

VC.Colors = 16

END IF

VC.Atribs = 4

VC.XPix = 319

VC.YPix = 199

VC.TCOL = 40

VC.TROW = 25

CASE 2 ' Two-color medium-res graphics for CGA, EGA, VGA,
and MCGA

IF BestMode = CGA OR BestMode = MCGA THEN

VC.Colors = 0

ELSE

VC.Colors = 16

END IF

VC.Atribs = 2

VC.XPix = 639

SUBSTITUTE SHEET

2184992

114

MenuSys.Bas - QBasic PopDown Menu System
SUB SetConfig

```
    VC.YPix = 199
    VC.TCOL = 80
    VC.TROW = 25
CASE 3    ' Two-color high-res graphics for Hercules
    VC.Colors = 0
    VC.Atribs = 2
    VC.XPix = 719
    VC.YPix = 347
    VC.TCOL = 80
    VC.TROW = 25
CASE 7    ' 16-color medium-res graphics for EGA and VGA
    VC.Colors = 16
    VC.Atribs = 16
    VC.XPix = 319
    VC.YPix = 199
    VC.TCOL = 40
    VC.TROW = 25
CASE 8    ' 16-color high-res graphics for EGA and VGA
    VC.Colors = 16
    VC.Atribs = 16
    VC.XPix = 639
    VC.YPix = 199
    VC.TCOL = 80
    VC.TROW = 25
CASE 9    ' 16- or 4-color very high-res graphics for EGA
and VGA
    VC.Colors = 64
    IF BestMode = EGA64 THEN VC.Atribs = 4 ELSE VC.Atribs
= 16
    VC.XPix = 639
    VC.YPix = 349
    VC.TCOL = 80
    VC.TROW = 25
```

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB SetConfig

CASE 10 ' Two-color high-res graphics for EGA or VGA
monochrome

VC.Colors = 0

VC.Atribs = 2

VC.XPix = 319

VC.YPix = 199

VC.TCOL = 80

VC.TROW = 25

CASE 11 ' Two-color very high-res graphics for VGA and
MCGA

' Note that for VGA screens 11, 12, and 13, more colors
are

' available, depending on how the colors are mixed.

VC.Colors = 216

VC.Atribs = 2

VC.XPix = 639

VC.YPix = 479

VC.TCOL = 80

VC.TROW = 30

CASE 12 ' 16-color very high-res graphics for VGA

VC.Colors = 216

VC.Atribs = 16

VC.XPix = 639

VC.YPix = 479

VC.TCOL = 80

VC.TROW = 30

CASE 13 ' 256-color medium-res graphics for VGA and MCGA

VC.Colors = 216

VC.Atribs = 256

VC.XPix = 639

VC.YPix = 479

VC.TCOL = 40

VC.TROW = 25

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB SetConfig

CASE ELSE

VC.Colors = 16

VC.Atribs = 16

VC.XPix = 0

VC.YPix = 0

VC.TCOL = 80

VC.TROW = 25

VC.Scrn = 0

EXIT SUB

END SELECT

VC.Scrn = mode

END SUB

FUNCTION StrTok\$ (Srce\$, Delim\$)

STATIC Start%, SaveStr\$

' If first call, make a copy of the string.

IF Srce\$ <> "" THEN

Start% = 1: SaveStr\$ = Srce\$

END IF

BegPos% = Start%: Ln% = LEN(SaveStr\$)

' Look for start of a token (character that isn't delimiter).

WHILE BegPos% <= Ln% AND INSTR(Delim\$, MID\$(SaveStr\$, BegPos%,
1)) <> 0

BegPos% = BegPos% + 1

WEND

' Test for token start found.

IF BegPos% > Ln% THEN

StrTok\$ = "": EXIT FUNCTION

END IF

' Find the end of the token.

EndPos% = BegPos%

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB SetConfig

```
    WHILE EndPos% <= Ln% AND INSTR(Delim$, MID$(SaveStr$, EndPos%,  
1)) = 0
```

```
        EndPos% = EndPos% + 1
```

```
    WEND
```

```
    StrTok$ = MID$(SaveStr$, BegPos%, EndPos% - BegPos%)
```

```
    ' Set starting point for search for next token.
```

```
    Start% = EndPos%
```

END FUNCTION

'\$Page \$SubTitle:'SUB TestGrModes'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB TestGrModes

SUB TestGrModes

 SHARED VC AS Config, InitRows AS INTEGER, BestMode AS INTEGER,
 Available AS STRING

 DIM i AS INTEGER, mode AS INTEGER, a\$

 'test all modes

 FOR i = 1 TO LEN(Available)

 a\$ = MID\$(Available, i, 1)

 IF a\$ >= "1" AND a\$ <= "9" THEN

 mode = VAL(a\$)

 ELSE

 mode = (ASC(a\$) - ASC("A") + 1) + 9

 END IF

 SetConfig mode

 SCREEN VC.Scrn

 PRINT "Vc.Scrn"; VC.Scrn

 PRINT "Vc.Colors"; VC.Colors

 PRINT "Vc.Atribs"; VC.Atribs

 PRINT "Vc.XPix"; VC.XPix

 PRINT "Vc.YPix"; VC.YPix

 PRINT "Vc.TCOL"; VC.TCOL

 PRINT "Vc.TROW"; VC.TROW

 PRINT "InitRows"; InitRows

 PRINT "BestMode"; BestMode

 PRINT "Available"; Available

 INPUT a\$

 SCREEN 0

 WIDTH 80, InitRows

 NEXT i

END SUB

'\$Page \$SubTitle:'SUB UnPop'

SUBSTITUTE SHEET

MenuSys.Bas - QBasic PopDown Menu System

SUB UnPop

SUB UnPop

 SHARED PopBuf() AS INTEGER, ForGnd AS INTEGER, BakGnd AS
INTEGER

 DIM r AS INTEGER, c AS INTEGER, i AS INTEGER

 i = 3

 FOR r = PopBuf(1, 1) TO PopBuf(2, 1)

 FOR c = PopBuf(1, 0) TO PopBuf(2, 0)

 LOCATE r, c

 COLOR (PopBuf(i, 1) AND &HF), (PopBuf(i, 1) AND &HF0) \ 16

 PRINT CHR\$(PopBuf(i, 0));

 i = i + 1

 NEXT c

 NEXT r

 COLOR ForGnd, BakGnd

END SUB

43949 Bytes Available

26833 Bytes Free

0 Warning Error(s)

0 Severe Error(s)

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

'\$Title:'VentLib.Bas - Ohio 7000 Ventilator Control'
\$LineSize:112'

'VentLib
'Oct 2/93
'Chris McLennan

'MenuSys.Bas
DECLARE FUNCTION Fuse\$ (x AS SINGLE, wide AS INTEGER, dec AS
INTEGER)
DECLARE FUNCTION LPad\$ (s\$, length AS INTEGER)
DECLARE FUNCTION StrTok\$ (Srce\$, Delim\$)

DECLARE SUB BrList (txt\$(), NumLin AS INTEGER)
DECLARE SUB MsgBox (txt\$(), NumLin AS INTEGER, Ky\$, RestoreScrn
AS INTEGER, Pause AS INTEGER)
DECLARE SUB SelectFile (FileName\$)
DECLARE SUB SetConfig (Mode AS INTEGER)
DECLARE SUB UnPop ()

'DasLib.Bas
DECLARE FUNCTION A2dToVolt! (A2dVal AS INTEGER)
DECLARE FUNCTION VoltToD2a\$ (Volt!)
DECLARE SUB DasMode0 ()
DECLARE SUB DasMode1 (ChLow AS INTEGER, ChHigh AS INTEGER)
DECLARE SUB DasMode2 (NextCh AS INTEGER, ChLow AS INTEGER,
ChHigh AS INTEGER)
DECLARE SUB DasMode3 (A2dData AS INTEGER, A2dCh AS INTEGER)
DECLARE SUB DasMode4 (NumPts AS INTEGER)
DECLARE SUB DasMode5 (NumPts AS INTEGER, Cycle AS INTEGER)
DECLARE SUB DasMode6 (NumPts AS INTEGER, Cycle AS INTEGER)
DECLARE SUB DasMode7 ()
DECLARE SUB DasMode8 (Op AS INTEGER, Status AS INTEGER, Count
AS INTEGER)
DECLARE SUB DasMode9 (NumPts AS INTEGER, StartPt AS INTEGER)

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

```

DECLARE SUB DasModel5 (D2aCh AS INTEGER, D2aData AS INTEGER)
DECLARE SUB DasModel6 (D2aDat0 AS INTEGER, D2aDat1 AS INTEGER)
DECLARE SUB DasModel7 (Rate!)

'VentLib.Bas
DECLARE SUB Demo ()
DECLARE SUB InitVentLib ()
DECLARE SUB LoopTest ()
DECLARE SUB LoopTest5 ()
DECLARE SUB LoopTest6 ()
DECLARE SUB OhioFn ()
DECLARE SUB TestD2a ()
DECLARE SUB TestMode5 ()
DECLARE SUB TestMode6 ()
DECLARE SUB VentLoop ()
DECLARE SUB VentPlotGrid ()
DECLARE FUNCTION AvgA2dChan% (Count AS INTEGER, Chan AS
INTEGER, Period AS INTEGER, NumCh AS
INTEGER, NumPts AS INTEGER)
DECLARE FUNCTION RateToV! (Opt AS INTEGER, Volt!)
DECLARE FUNCTION VoltToV! (Opt AS INTEGER, Volt!)

'$INCLUDE: 'DasLib.Bil'
'----- DasLib.Bil -----

'Das16 Common Area

'Das16 Parameters MUST be in common

CONST A2dBuffSize = 2000

DIM Dio(4) AS INTEGER
DIM A2dBuff(2000) AS INTEGER
DIM A2dChNum(2000) AS INTEGER

```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

```
COMMON Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER
COMMON BiPolar AS INTEGER, DasChan AS INTEGER, VoltFS AS SINGLE
COMMON A2dBuff() AS INTEGER           'Buffer for DasMode4
COMMON A2dChNum() AS INTEGER          'Buffer for A/D channel
info
```

```
'for Mode 5 & 6, Quick Basic manages the Data Segment
```

```
'$DYNAMIC
```

```
'Dynamic Common is dimensioned AFTER the COMMON specification
```

```
CONST A2dBuf2Size = 2000
```

```
COMMON A2dBuf2() AS INTEGER           'Buffer for A/D data
```

```
'$STATIC
```

```
'----- End DasLib.Bil -----
```

```
'$INCLUDE: 'MenuSys.Bil'
```

```
'----- MenuSys.Bil-----
```

```
' Constants for best available screen mode
```

```
CONST VGA = 12
```

```
CONST MCGA = 13
```

```
CONST EGA256 = 9
```

```
CONST EGA64 = 8
```

```
CONST MONO = 10
```

```
CONST HERC = 3
```

```
CONST CGA = 1
```

```
' User-defined type to hold information about the mode
```

```
TYPE Config
```

```
    Scrn      AS INTEGER
```

```
    Colors    AS INTEGER
```

```
    Atribs    AS INTEGER
```

```
    XPix      AS INTEGER
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

```

    YPix      AS INTEGER
    TCOL      AS INTEGER
    TROW      AS INTEGER
END TYPE

```

```

'*****Graphics Modes

```

```

COMMON VC AS Config, InitRows AS INTEGER, BestMode AS INTEGER,
Available AS STRING

```

```

'****Menu system

```

```

'$DYNAMIC

```

```

COMMON PopBuf() AS INTEGER, CurPath$, ForGnd AS INTEGER, Bakgnd
AS INTEGER

```

```

'$STATIC

```

```

'----- End MenuSys.Bil-----

```

```

'$INCLUDE: 'VentLib.Bil'

```

```

'----- VentLib.Bil -----

```

```

'VentLib Dynamic Common declaration

```

```

TYPE VentDat

```

```

'Structure to hold modulation data

```

```

    Time AS SINGLE

```

```

    Mode AS STRING * 1

```

```

    Modulation AS SINGLE

```

```

END TYPE

```

```

CONST VentArrSize = 2000

```

```

COMMON VentArrNum AS INTEGER
VentArr()

```

```

'# of items loaded in

```

```

COMMON VentDatFn$

```

```

'Data File loaded

```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

'\$DYNAMIC

'Dynamic Common is dimensioned AFTER the COMMON specification

COMMON VentArr() AS VentDat 'Buffer for Ventilator
modulation data

'\$STATIC

'----- End VentLib.Bi1 -----

'\$INCLUDE: 'DasLib.Bi2'

'----- DasLib.Bi2 -----

DIM A2dBuf2(A2dBuf2Size) AS INTEGER

'----- End DasLib.Bi2 -----

'\$INCLUDE: 'MenuSys.Bi2'

'----- MenuSys.Bi2-----

DIM PopBuf(2002, 1) AS INTEGER 'PopUp Buffer

'----- End MenuSys.Bi2-----

'\$INCLUDE: 'VentLib.Bi2'

'----- VentLib.Bi2 -----

DIM VentArr(VentArrSize) AS VentDat 'Buffer for Ventilator
modulation data

'----- End VentLib.Bi2 -----

'\$Page \$SubTitle: 'FUNCTION AvgA2dChan%

SUBSTITUTE SHEET

2184992

125

VentLib.Bas - Ohio 7000 Ventilator Control

FUNCTION AvgA2dChan%

FUNCTION AvgA2dChan% (Count AS INTEGER, Chan AS INTEGER, Period
AS INTEGER, NumCh AS INTEGER, N
umPts AS INTEGER)

'Function returns the average of an A/D channel over the
previous Period

'of observations, ending at position Count.

'Data is input into A2dBuff() and A2dChNum() before call by
DasMode9.

'Count is provided by DasMode8, before call to DasMode9.

' is in range [0..(NumCh*NumPts)-1].

'Chan is the channel # to average.

'Period is the # of previous observations to average over.

'NumCh channels of data is in data arrays.

'NumPts is the number of points per channel.

'The index pointers may wrap around.

SHARED A2dBuff() AS INTEGER	'A/D data buffer
SHARED A2dChNum() AS INTEGER	'channel # buffer

DIM ptr AS INTEGER, i AS INTEGER
DIM sum

IF Count < 0 OR Count > (NumCh * NumPts) - 1 THEN
PRINT "Count out of range in AvgA2dChan"
DasMode7
STOP
END IF

ptr = Count - 1 'use last conversion as end

'scan backwards until desired channel is located
DO

SUBSTITUTE SHEET

126 2184992

VentLib.Bas - Ohio 7000 Ventilator Control

FUNCTION AvgA2dChan%

```
    IF ptr < 0 THEN ptr = ptr + (NumCh * NumPts)      'protect
wrap arounds
```

```
    IF A2dChNum(ptr) = Chan THEN EXIT DO
    ptr = ptr - 1
LOOP
```

```
sum = 0
```

```
FOR i = 1 TO Period
    IF ptr < 0 THEN ptr = ptr + (NumCh * NumPts)      'protect
wrap arounds
```

```
    IF A2dChNum(ptr) = Chan THEN
        sum = sum + A2dBuff(ptr)
```

```
    END IF
    ptr = ptr - NumCh
NEXT i
```

```
AvgA2dChan% = sum / Period
```

```
END FUNCTION
```

```
'$Page $SubTitle:'SUB BrPrnData'
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB BrPrnData

SUB BrPrnData

'Browse Ventilator modulation data

SHARED CurPath\$, ForGnd AS INTEGER, Bakgnd AS INTEGER

SHARED VentArrNum AS INTEGER '# of items loaded in
VentArr()

SHARED VentDatFn\$ 'Data File loaded

SHARED VentArr() AS VentDat 'Buffer for Ventilator
modulation data

DIM i AS INTEGER, Ky\$

IF VentArrNum = 0 THEN

REDIM Msg\$(3)

Msg\$(0) = "Browse Data"

Msg\$(1) = ""

Msg\$(2) = "NO data loaded"

Msg\$(3) = ""

MsgBox Msg\$(), 3, Ky\$, -1, -1 'unpop, pause

ERASE Msg\$

EXIT SUB 'nothing to do

END IF

'entertain the Human

REDIM Msg\$(3)

Msg\$(0) = "Browse Data"

Msg\$(1) = ""

Msg\$(2) = "Formating List..."

Msg\$(3) = ""

SUBSTITUTE SHEET

2184992

128

VentLib.Bas - Ohio 7000 Ventilator Control
SUB BrPrnData

MsgBox Msg\$(), 3, Ky\$, 0, 0 'no unpop, no pause
ERASE Msg\$

'Create list to browse
REDIM txt\$(VentArrNum)

'Generate header
txt\$(0) = "File: " + VentDatFn\$
i = 80 - LEN(txt\$(0))
txt\$(0) = txt\$(0) + LPad\$("ESC - Exit", i)

FOR i = 1 TO VentArrNum
 txt\$(i) = LPad\$(STR\$(i), 5)
 txt\$(i) = txt\$(i) + FUse\$(VentArr(i).Time, 10, 2)
 txt\$(i) = txt\$(i) + LPad\$(LEFT\$(VentArr(i).Mode, 1), 5)
 txt\$(i) = txt\$(i) + FUse\$(VentArr(i).Modulation, 10, 2)
NEXT i

UnPop 'Clear Message box

BrList txt\$(), VentArrNum

ERASE txt\$

END SUB

'\$Page \$SubTitle:'SUB demo'

SUBSTITUTE SHEET

SUB demo

```

SHARED Dio() AS INTEGER, Md AS INTEGER, Flag AS INTEGER
SHARED BiPolar AS INTEGER, DasChan AS INTEGER, VoltFS AS SINGLE
SHARED A2dBuff() AS INTEGER 'Buffer for DasMode4
SHARED A2dBuf2() AS INTEGER

```

SUBSTITUTE SHEET

130

VentLib.Bas - Ohio 7000 Ventilator Control
SUB demo

TestMode5

INPUT "mode 6 test"; a\$
TestMode6

TestD2a

END SUB

'\$Page \$SubTitle:'SUB InitDas16'

SUBSTITUTE SHEET

2184992

131

VentLib.Bas - Ohio 7000 Ventilator Control
SUB InitDas16

SUB InitDas16

'General Reset

DasMode0

'Init Driver

DasMode7

'Clear DMA or Interrupt operation

DasMode16 0, 0

'D/A's to zero

END SUB

'\$Page \$SubTitle:'SUB InitVentLib'

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB InitVentLib

SUB InitVentLib

 SHARED VentArrNum AS INTEGER '# of items loaded in
VentArr()
 SHARED VentDatFn\$ 'Data File loaded
 SHARED VentArr() AS VentDat 'Buffer for Ventilator
modulation data

VentArrNum = 0
VentDatFn\$ = ""

END SUB

'\$Page \$SubTitle:'SUB LoadVentDatPrn'

SUBSTITUTE SHEET

133

VentLib.Bas - Ohio 7000 Ventilator Control

SUB LoadVentDatPrn

SUB LoadVentDatPrn

'Load a .PRN format file created by a spreadsheet to emulate
real data'

'	Time	Rate/Volume	Modulation
'	0	R	6
'	0	V	2
'	4.23569	R	36.81015
'	13.60336	R	13.74554
'	20.26031	V	28.91379
'	25.65653	V	14.68135

SHARED CurPath\$, ForGnd AS INTEGER, Bakgnd AS INTEGER

SHARED VentArrNum AS INTEGER '# of items loaded in
VentArr()

SHARED VentDatFn\$ 'Data File loaded

SHARED VentArr() AS VentDat 'Buffer for Ventilator
modulation data

DIM Ky\$, p\$

DIM fp AS INTEGER, i AS INTEGER

REDIM Msg\$(5)

Msg\$(0) = "Import VENTILATOR data"

Msg\$(1) = ""

Msg\$(2) = "Select an ASCII data file to load"

Msg\$(3) = "as a Ventilator control data file."

Msg\$(4) = ""

Msg\$(5) = "Press any key to continue"

MsgBox Msg\$(), 5, Ky\$, -1, -1 'unpop, pause

ERASE Msg\$

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoadVentDatPrn

IF Ky\$ = CHR\$(27) THEN EXIT SUB 'ESC pressed

'Try to open a file

SelectFile VentDatFn\$

IF VentDatFn\$ = "" THEN

'nothing selected

InitVentLib

'clear common buffer

EXIT SUB

END IF

'Display status box while loading

REDIM Msg\$(3)

Msg\$(0) = "Import VENTILATOR data"

Msg\$(1) = ""

Msg\$(2) = "Loading: " + CurPath\$ + VentDatFn\$

Msg\$(3) = ""

MsgBox Msg\$(), 3, Ky\$, 0, 0 'no unpop, no pause

ERASE Msg\$

fp = FREEFILE

OPEN CurPath\$ + VentDatFn\$ FOR INPUT AS #fp

VentArrNum = 0

'init # loaded

i = 1

DO WHILE NOT EOF(1)

LINE INPUT #1, p\$

p\$ = UCASE\$(LTRIM\$(RTRIM\$(p\$)))

IF LEN(p\$) > 0 THEN

'Parse input line

VentArr(0).Time = VAL(StrTok\$(p\$, " "))

SUBSTITUTE SHEET

2184992

135

VentLib.Bas - Ohio 7000 Ventilator Control

SUB LoadVentDatPrn

```
VentArr(0).Mode = LEFT$(StrTok$("", " "), 1)
VentArr(0).Modulation = VAL(StrTok$("", " "))
```

```
'test if valid - see format above
```

```
IF VentArr(0).Time >= 0 AND (VentArr(0).Mode = "R" OR
VentArr(0).Mode = "V") AND VentArr(
0).Modulation > 0 THEN
```

```
VentArr(i).Time = VentArr(0).Time
VentArr(i).Mode = VentArr(0).Mode
VentArr(i).Modulation = VentArr(0).Modulation
VentArrNum = i
i = i + 1
```

```
END IF
```

```
'test if buffer exceeded
```

```
IF i > VentArrSize THEN EXIT DO
END IF
```

```
LOOP
```

```
CLOSE #fp
```

```
UnPop
```

```
'Clear "Loading" Message - no unpop
```

```
IF VentArrNum = 0 THEN
```

```
VentDatFn$ = ""
```

```
EXIT SUB
```

```
END IF
```

```
REDIM Msg$(5)
```

```
Msg$(0) = "Import VENTILATOR data"
```

```
Msg$(1) = ""
```

```
Msg$(2) = "File: " + CurPath$ + VentDatFn$
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoadVentDatPrn

Msg\$(3) = STR\$(VentArrNum) + " obs. loaded"

Msg\$(4) = ""

Msg\$(5) = "Press any key to continue"

MsgBox Msg\$(), 5, Ky\$, -1, -1 'unpop, pause

ERASE Msg\$

END SUB

'\$Page \$SubTitle:'SUB LoopTest'

SUBSTITUTE SHEET

2184992

137

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest

SUB LoopTest

'Closed Loop Test of Ohio Ventilator
'R decreases modulated rate, T increases
'V decreases modulated volume, B increases
'Uses Mode4 Das16 operation

SHARED A2dBuff() AS INTEGER 'Buffer for DasMode4

DIM cmd\$, Avg(3), i AS INTEGER, VolCtrl, RateCtrl, VolMod,
RateMod

DIM VolOut, RateOut, VolD2a, RateD2a

CLS

PRINT "ESC to terminate Closed Loop Test"
PRINT "R/T = -/+ Rate V/B = -/+ Volume"

LOCATE 5, 15: PRINT "---Volume--"
LOCATE 6, 15: PRINT "L/Min Volt";

LOCATE 5, 30: PRINT "----Rate----"
LOCATE 6, 30: PRINT "Bth/M Volt";

LOCATE 7, 1: PRINT "Baseline";
LOCATE 8, 1: PRINT "Modulation";
LOCATE 9, 1: PRINT "D/A Drive";
LOCATE 10, 1: PRINT "Output";

DasModel 0, 3 'Setup Das16 for channels 0-3
DasModel7 (100) '100 Hz/Ch sampling rate

'Initialize keyboard volume and rate modulation values
VolMod = 2
RateMod = 6

SUBSTITUTE SHEET

138

2184992

VentLib.Bas - Ohio 7000 Ventilator Control

SUB LoopTest

DO

```

      'Use Mode4 to get 100 mS of Data - speed medium
      DasMode4 (4 * 10)          '10 pts per channel @ 10 mS/Pt =
100 mS

```

```

      'Cal avg of the volume and rate control settings

```

```

      Avg(0) = 0          'Volume

```

```

      Avg(2) = 0          'Rate

```

```

      FOR i = 0 TO (10 * 4) - 1 STEP 4

```

```

          Avg(0) = Avg(0) + A2dBuff(i)

```

```

          Avg(2) = Avg(2) + A2dBuff(i + 2)

```

```

      NEXT i

```

```

      Avg(0) = A2dToVolt(CINT(Avg(0) / 10!))

```

```

      VolCtrl = VolToV(1, Avg(0))

```

```

      Avg(2) = A2dToVolt(CINT(Avg(2) / 10!))

```

```

      RateCtrl = RateToV(1, Avg(2))

```

```

      'Report Control settings to screen

```

```

      LOCATE 7, 15

```

```

      PRINT USING "##.## #.###"; VolCtrl; Avg(0);

```

```

      LOCATE 7, 30

```

```

      PRINT USING "##.## #.###"; RateCtrl; Avg(2);

```

```

      cmd$ = INKEY$

```

```

      SELECT CASE UCASE$(cmd$)

```

```

          CASE "R"

```

```

              IF RateMod > 6 THEN RateMod = RateMod - 1

```

```

          CASE "T"

```

```

              IF RateMod < 40 THEN RateMod = RateMod + 1

```

SUBSTITUTE SHEET

2184992

139

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest

```
CASE "V"
    IF VolMod > 2 THEN VolMod = VolMod - 1
CASE "B"
    IF VolMod < 30 THEN VolMod = VolMod + 1

CASE CHR$(27)
    EXIT DO
CASE ""
CASE ELSE
    BEEP
END SELECT

'Display the Volume and Rate Modulations
LOCATE 8, 15
PRINT USING "##.## #.###"; VolMod; VolToV(2, VolMod);
LOCATE 8, 30
PRINT USING "##.## #.###"; RateMod; RateToV(2, RateMod);

'if the Modulation level is greater than the set point,
'increase the drive to the ventilator

VolD2a = 0
IF VolMod > VolCtrl THEN
    VolD2a = VolToV(2, VolMod) - VolToV(2, VolCtrl)
END IF

RateD2a = 0
IF RateMod > RateCtrl THEN
    RateD2a = RateToV(2, RateMod) - RateToV(2, RateCtrl)
END IF

'Output modulation to D/A
DasModel6 VoltToD2a%(VolD2a), VoltToD2a%(RateD2a)
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest

```
'Display the Volume and Rate D/A Modulations levels
LOCATE 9, 15
PRINT USING "      ###"; VolD2a;
LOCATE 9, 30
PRINT USING "      ###"; RateD2a;

'Use Mode4 to get 100 mS of Data - speed medium
DasMode4 (4 * 10)          '10 pts per channel @ 10 mS/Pt =
100 mS

'Cal avg of the volume output and rate output levels
Avg(1) = 0                'Volume Output
Avg(3) = 0                'Rate Output
FOR i = 0 TO (10 * 4) - 1 STEP 4
    Avg(1) = Avg(1) + A2dBuff(i + 1)
    Avg(3) = Avg(3) + A2dBuff(i + 3)
NEXT i

Avg(1) = A2dToVolt(CINT(Avg(1) / 10!))
VolOut = VolToV(1, Avg(1))
Avg(3) = A2dToVolt(CINT(Avg(3) / 10!))
RateOut = RateToV(1, Avg(3))
'Report Outputs to screen
LOCATE 10, 15
PRINT USING "##.## ###"; VolOut; Avg(1);
LOCATE 10, 30
PRINT USING "##.## ###"; RateOut; Avg(3);
LOOP UNTIL cmd$ = CHR$(27)

'Zero D/A's
DasModel6 0, 0
END SUB

'$Page $SubTitle:'SUB LoopTest5'
```

SUBSTITUTE SHEET

2184992

141

VentLib.Bas - Ohio 7000 Ventilator Control

SUB LoopTest5

SUB LoopTest5

'Uses Mode5 Das16 operation

'Closed Loop Test of Ohio Ventilator

'R decreases modulated rate, T increases

'V decreases modulated volume, B increases

SHARED A2dBuff() AS INTEGER 'Buffer for DasMode4

DIM cmd\$, Avg(3), i AS INTEGER, VolCtrl, RateCtrl, VolMod,
RateMod

DIM VolOut, RateOut, VolD2a, RateD2a

DIM Op AS INTEGER, Status AS INTEGER, Count AS INTEGER

CLS

PRINT "ESC to terminate Closed Loop Test"

PRINT "R/T = -/+ Rate V/B = -/+ Volume"

LOCATE 5, 15: PRINT "----Volume--"

LOCATE 6, 15: PRINT "L/Min Volt";

LOCATE 5, 30: PRINT "----Rate----

LOCATE 6, 30: PRINT "Bth/M Volt";

LOCATE 7, 1: PRINT "Baseline";

LOCATE 8, 1: PRINT "Modulation";

LOCATE 9, 1: PRINT "D/A Drive";

LOCATE 10, 1: PRINT "Output";

DasMode1 0, 3 'Setup Das16 for channels 0-3

DasMode17 (100) '100 Hz/Ch sampling rate

'Initialize keyboard volume and rate modulation values

SUBSTITUTE SHEET

2184992

142

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest5

VolMod = 2
RateMod = 6

DO

'Use Mode5 to get 100 mS of Data - speed medium - Background
DasMode5 4 * 10, 0 'get 10 pts per channel, only one
cycle

'watch Count until done
DO

DasMode8 Op, Status, Count
LOOP UNTIL Status = 0

DasMode9 10 * 4, 0 'Xfer to A2dBuff()

'Cal avg of the volume and rate control settings
Avg(0) = 0 'Volume
Avg(2) = 0 'Rate
FOR i = 0 TO (10 * 4) - 1 STEP 4
Avg(0) = Avg(0) + A2dBuff(i)
Avg(2) = Avg(2) + A2dBuff(i + 2)
NEXT i

Avg(0) = A2dToVolt(CINT(Avg(0) / 10!))
VolCtrl = VolToV(1, Avg(0))
Avg(2) = A2dToVolt(CINT(Avg(2) / 10!))
RateCtrl = RateToV(1, Avg(2))

'Report Control settings to screen
LOCATE 7, 15
PRINT USING "##.## #.###"; VolCtrl; Avg(0);
LOCATE 7, 30

SUBSTITUTE SHEET

2184992

143

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest5

PRINT USING "##.## #.###"; RateCtrl; Avg(2);

cmd\$ = INKEY\$

SELECT CASE UCASE\$(cmd\$)

CASE "R"

IF RateMod > 6 THEN RateMod = RateMod - 1

CASE "T"

IF RateMod < 40 THEN RateMod = RateMod + 1

CASE "V"

IF VolMod > 2 THEN VolMod = VolMod - 1

CASE "B"

IF VolMod < 30 THEN VolMod = VolMod + 1

CASE CHR\$(27)

EXIT DO

CASE ""

CASE ELSE

BEEP

END SELECT

'Display the Volume and Rate Modulations

LOCATE 8, 15

PRINT USING "##.## #.###"; VolMod; VolToV(2, VolMod);

LOCATE 8, 30

PRINT USING "##.## #.###"; RateMod; RateToV(2, RateMod);

'if the Modulation level is greater than the set point,
'increase the drive to the ventilator

VolD2a = 0

IF VolMod > VolCtrl THEN

VolD2a = VolToV(2, VolMod) - VolToV(2, VolCtrl)

END IF

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest5

```
RateD2a = 0
IF RateMod > RateCtrl THEN
    RateD2a = RateToV(2, RateMod) - RateToV(2, RateCtrl)
END IF
```

```
'Output modulation to D/A
DasMode16 VoltToD2a%(VolD2a), VoltToD2a%(RateD2a)
```

```
'Display the Volume and Rate D/A Modulations levels
LOCATE 9, 15
PRINT USING "      ###"; VolD2a;
LOCATE 9, 30
PRINT USING "      ###"; RateD2a;
```

```
'Use Mode5 to get 100 mS of Data - speed medium - Background
DasMode5 4 * 10, 0      'get 10 pts per channel, only one
cycle
```

```
'watch Count until done
DO
    DasMode8 Op, Status, Count
LOOP UNTIL Status = 0
```

```
DasMode9 10 * 4, 0      'Xfer to A2dBuff()
```

```
'Cal avg of the volume output and rate output levels
Avg(1) = 0      'Volume Output
Avg(3) = 0      'Rate Output
FOR i = 0 TO (10 * 4) - 1 STEP 4
    Avg(1) = Avg(1) + A2dBuff(i + 1)
    Avg(3) = Avg(3) + A2dBuff(i + 3)
NEXT i
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB LoopTest5

Avg(1) = A2dToVolt(CINT(Avg(1) / 10!))

VolOut = VolToV(1, Avg(1))

Avg(3) = A2dToVolt(CINT(Avg(3) / 10!))

RateOut = RateToV(1, Avg(3))

'Report Outputs to screen

LOCATE 10, 15

PRINT USING "##.## #.###"; VolOut; Avg(1);

LOCATE 10, 30

PRINT USING "##.## #.###"; RateOut; Avg(3);

LOOP UNTIL cmd\$ = CHR\$(27)

DasModel6 0, 0

'Zero D/A's

DasMode7

'Disable Interrupt

END SUB

'\$Page \$SubTitle:'SUB LoopTest6'

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB LoopTest6

SUB LoopTest6

'Uses Mode6 Das16 operation

'Closed Loop Test of Ohio Ventilator

'R decreases modulated rate, T increases

'V decreases modulated volume, B increases

SHARED A2dBuff() AS INTEGER 'Buffer for DasMode4

DIM cmd\$, Avg(3), i AS INTEGER, VolCtrl, RateCtrl, VolMod,
RateMod

DIM VolOut, RateOut, VolD2a, RateD2a

DIM Op AS INTEGER, Status AS INTEGER, Count AS INTEGER

CLS

PRINT "ESC to terminate Closed Loop Test"

PRINT "R/T = -/+ Rate V/B = -/+ Volume"

LOCATE 5, 15: PRINT "---Volume--"

LOCATE 6, 15: PRINT "L/Min Volt";

LOCATE 5, 30: PRINT "----Rate---

LOCATE 6, 30: PRINT "Bth/M Volt";

LOCATE 7, 1: PRINT "Baseline";

LOCATE 8, 1: PRINT "Modulation";

LOCATE 9, 1: PRINT "D/A Drive";

LOCATE 10, 1: PRINT "Output";

DasModel 0, 3 'Setup Das16 for channels 0-3

DasModel7 (100) '100 Hz/Ch sampling rate

'Initialize keyboard volume and rate modulation values

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest6

VolMod = 2

RateMod = 6

DO

'Use Mode6 to get 100 mS of Data - speed high - Background
DasMode6 4 * 10, 0 'get 10 pts per channel, only one
cycle

'watch Count until done

DO

DasMode8 Op, Status, Count
LOOP UNTIL Status = 0

DasMode9 10 * 4, 0 'Xfer to A2dBuff()

'Cal avg of the volume and rate control settings

Avg(0) = 0 'Volume

Avg(2) = 0 'Rate

FOR i = 0 TO (10 * 4) - 1 STEP 4

Avg(0) = Avg(0) + A2dBuff(i)

Avg(2) = Avg(2) + A2dBuff(i + 2)

NEXT i

Avg(0) = A2dToVolt(CINT(Avg(0) / 10!))

VolCtrl = VolToV(1, Avg(0))

Avg(2) = A2dToVolt(CINT(Avg(2) / 10!))

RateCtrl = RateToV(1, Avg(2))

'Report Control settings to screen

LOCATE 7, 15

PRINT USING "##.## #.###"; VolCtrl; Avg(0);

LOCATE 7, 30

PRINT USING "##.## #.###"; RateCtrl; Avg(2);

SUBSTITUTE SHEET

2184992

148

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest6

cmd\$ = INKEY\$

SELECT CASE UCASE\$(cmd\$)

CASE "R"

IF RateMod > 6 THEN RateMod = RateMod - 1

CASE "T"

IF RateMod < 40 THEN RateMod = RateMod + 1

CASE "V"

IF VolMod > 2 THEN VolMod = VolMod - 1

CASE "B"

IF VolMod < 30 THEN VolMod = VolMod + 1

CASE CHR\$(27)

EXIT DO

CASE ""

CASE ELSE

BEEP

END SELECT

'Display the Volume and Rate Modulations

LOCATE 8, 15

PRINT USING "##.## #.###"; VolMod; VolToV(2, VolMod);

LOCATE 8, 30

PRINT USING "##.## #.###"; RateMod; RateToV(2, RateMod);

'if the Modulation level is greater than the set point,

'increase the drive to the ventilator

VolD2a = 0

IF VolMod > VolCtrl THEN

VolD2a = VolToV(2, VolMod) - VolToV(2, VolCtrl)

END IF

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest6

```
RateD2a = 0
IF RateMod > RateCtrl THEN
    RateD2a = RateToV(2, RateMod) - RateToV(2, RateCtrl)
END IF

'Output modulation to D/A
DasMode16 VoltToD2a$(VolD2a), VoltToD2a$(RateD2a)

'Display the Volume and Rate D/A Modulations levels
LOCATE 9, 15
PRINT USING "      ###"; VolD2a;
LOCATE 9, 30
PRINT USING "      ###"; RateD2a;

'Use Mode6 to get 100 mS of Data - speed high - Background
DasMode6 4 * 10, 0      'get 10 pts per channel, only one
cycle

'watch Count until done
DO
    DasMode8 Op, Status, Count
LOOP UNTIL Status = 0

DasMode9 10 * 4, 0      'Xfer to A2dBuff()

'Cal avg of the volume output and rate output levels
Avg(1) = 0      'Volume Output
Avg(3) = 0      'Rate Output
FOR i = 0 TO (10 * 4) - 1 STEP 4
    Avg(1) = Avg(1) + A2dBuff(i + 1)
    Avg(3) = Avg(3) + A2dBuff(i + 3)
NEXT i
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB LoopTest6

```
Avg(1) = A2dToVolt(CINT(Avg(1) / 10!))  
VolOut = VolToV(1, Avg(1))  
Avg(3) = A2dToVolt(CINT(Avg(3) / 10!))  
RateOut = RateToV(1, Avg(3))
```

```
'Report Outputs to screen  
LOCATE 10, 15  
PRINT USING "##.## #.###"; VolOut; Avg(1);  
LOCATE 10, 30  
PRINT USING "##.## #.###"; RateOut; Avg(3);
```

```
LOOP UNTIL cmd$ = CHR$(27)
```

```
DasModel6 0, 0          'Zero D/A's  
DasMode7               'Disable Interrupt
```

END SUB

'\$Page \$SubTitle:'SUB OhioFn'

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB OhioFn

SUB OhioFn

'VentFn4.Bas

'Sept 29/93

DIM i AS INTEGER, Volt, Vol, Rate

DasMode0

DasMode1 0, 3

DasMode7

CLS

LOCATE 4, 1

PRINT "Volt";

LOCATE 5, 1

PRINT "Calc";

LOCATE 3, 10

PRINT "-----Litre/Min-----";

LOCATE 3, 40

PRINT "-----Breath/Min-----";

DO

FOR ch = 0 TO 3

DasMode3 A2dData%, A2dCh%

' PRINT A2dCh%;

IF ch < 2 THEN

LOCATE 6, 10 + (ch * 15)

PRINT USING "#####"; A2dData%;

LOCATE 4, 10 + (ch * 15)

PRINT USING "#.###"; A2dToVolt(A2dData%);

IF ch = 1 THEN DasMode15 0,

VoltToD2a%(A2dToVolt(A2dData%))

SUBSTITUTE SHEET

152

VentLib.Bas - Ohio 7000 Ventilator Control
SUB OhioFn

LOCATE 5, 10 + (ch * 15)
PRINT USING "##.##"; VolToV(1, A2dToVolt(A2dData%));

ELSE

LOCATE 4, 10 + (ch * 15)
PRINT USING "#.###"; A2dToVolt(A2dData%);
LOCATE 5, 10 + (ch * 15)
PRINT USING "##.##"; RateToV(1, A2dToVolt(A2dData%));

END IF

NEXT ch

PRINT

a\$ = INKEY\$

LOOP UNTIL a\$ <> ""

EXIT SUB

FOR i = 2 TO 30

Volt = VolToV(2, CSNG(i))

Vol = VolToV(1, Volt)

PRINT USING "## #.### ####.#"; i; Volt; Vol

IF i MOD 10 = 0 THEN INPUT a\$

NEXT i

FOR i = 6 TO 40

Volt = RateToV(2, CSNG(i))

Rate = RateToV(1, Volt)

PRINT USING "## #.### ####.#"; i; Volt; Rate

IF i MOD 10 = 0 THEN INPUT a\$

NEXT i

END SUB

'\$Page \$SubTitle:'SUB PlotVentDat'

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB PlotVentDat

SUB PlotVentDat

 SHARED VentArrNum AS INTEGER '# of items loaded in
VentArr()
 SHARED VentDatFn\$ 'Data File loaded
 SHARED VentArr() AS VentDat 'Buffer for Ventilator
modulation data

DIM i AS INTEGER, init AS INTEGER
DIM oldY AS SINGLE

IF VentArrNum = 0 THEN

 REDIM Msg\$(3)
 Msg\$(0) = "Plot Data"
 Msg\$(1) = ""
 Msg\$(2) = "NO data loaded"
 Msg\$(3) = ""

 MsgBox Msg\$(), 3, Ky\$, -1, -1 'unpop, pause

 ERASE Msg\$

 EXIT SUB 'nothing to do

END IF

VentPlotGrid 'draw axis and define plot window

'Plot Volume

WINDOW (0, 0)-(VentArr(VentArrNum).Time, 30)

init = -1

FOR i = 1 TO VentArrNum

 IF VentArr(i).Mode = "V" THEN 'Volume data

 IF init THEN

 init = 0

SUBSTITUTE SHEET

2184992

154

VentLib.Bas - Ohio 7000 Ventilator Control

SUB PlotVentDat

```
        PSET (VentArr(i).Time, VentArr(i).Modulation), 1
        oldY = VentArr(i).Modulation
    ELSE
        LINE -(VentArr(i).Time, oldY), 1
        LINE -(VentArr(i).Time, VentArr(i).Modulation), 1
        oldY = VentArr(i).Modulation
    END IF
END IF
NEXT i

'Plot Rate
WINDOW (0, 0)-(VentArr(VentArrNum).Time, 40)
init = -1
FOR i = 1 TO VentArrNum
    IF VentArr(i).Mode = "R" THEN          'Rate data
        IF init THEN
            init = 0
            PSET (VentArr(i).Time, VentArr(i).Modulation), 2
            oldY = VentArr(i).Modulation
        ELSE
            LINE -(VentArr(i).Time, oldY), 2
            LINE -(VentArr(i).Time, VentArr(i).Modulation), 2
            oldY = VentArr(i).Modulation
        END IF
    END IF
NEXT i
DO
    a$ = INKEY$
LOOP UNTIL a$ <> ""
SCREEN 0, 0, 0
WIDTH 80
END SUB
```

'\$Page \$SubTitle:'FUNCTION RateToV'

SUBSTITUTE SHEET

155

2184992

VentLib.Bas - Ohio 7000 Ventilator Control

FUNCTION RateToV

FUNCTION RateToV (Opt AS INTEGER, Volt)

'Sept 29/93

'-----

'Function Library for Ohio 7000 Ventilator Control

'-----

'Based on observations of Breath/Min Rate control

' Y = a + bX

' mV = 18.06555 + 50.48908 * Rate R^2=0.999555

' V = 1.80655462E-02 + 5.04890756E-02 * Rate

R^2=0.999555

' Rate Regression Output:

' Constant 1.80655462E-02

' Std Err of Y Est 1.10777101E-02

' R Squared 9.99555207E-01

' No. of Observations 35

' Degrees of Freedom 33

'

' X Coefficient(s) 5.04890756E-02

' Std Err of Coef. 1.85402629E-04

DIM V, Rate, a, b

' V = 1.80655462E-02 + 5.04890756E-02 * Rate

R^2=0.999555

a = .0180655462#

b = .0504890756#

SELECT CASE Opt%

CASE 1 'V to Rate

V = Volt

SUBSTITUTE SHEET

156

2184992

VentLib.Bas - Ohio 7000 Ventilator Control

FUNCTION RateToV

 $\text{Rate} = (V - a) / b$ ' $x = (Y - a) / b$

RateToV = Rate

CASE 2 'Rate to V

Rate = Volt

 $V = a + b * \text{Rate}$ ' $y = a + bx$

RateToV = V

END SELECT

END FUNCTION

'\$Page \$SubTitle:'SUB TestD2a'

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB TestD2a

SUB TestD2a

DIM i AS INTEGER, a\$

INPUT "Channel 0 ramp"; a\$

FOR i = 0 TO 4095 STEP 2

DasModel5 0, i

NEXT i

DasModel5 0, 0

INPUT "Channel 1 ramp"; a\$

FOR i = 0 TO 4095 STEP 2

DasModel5 1, i

NEXT i

DasModel5 1, 0

INPUT "Dual Channel ramp"; a\$

FOR i = 0 TO 4095 STEP 2

DasModel6 i, i

NEXT i

DasModel6 0, 0

END SUB

'\$Page \$SubTitle:'SUB TestMode5'

SUBSTITUTE SHEET

2184992

158

VentLib.Bas - Ohio 7000 Ventilator Control
SUB TestMode5

SUB TestMode5

SHARED A2dBuff() AS INTEGER 'A/D data buffer
SHARED A2dChNum() AS INTEGER 'channel # buffer

DIM Op AS INTEGER, Status AS INTEGER, Count AS INTEGER, i AS
INTEGER

DasMode0 'init
DasMode1 0, 3 'scan 0 to 3
DasMode17 (100) 'sample @ 100 Hz
DasMode5 500, 0 'get 500 pts, only one cycle

DO

 DasMode8 Op, Status, Count
 PRINT Op, Status, Count
LOOP UNTIL Status = 0
PRINT Op, Status, Count

DasMode9 500, 0 'Xfer to A2dBuff(), A2dChNum()

FOR i = 0 TO 499 STEP 4
 PRINT A2dChNum(i); A2dBuff(i), A2dChNum(i + 1); A2dBuff(i +
1), A2dChNum(i + 2); A2dBuff(i
+ 2), A2dChNum(i + 3); A2dBuff(i + 3)
NEXT i

DasMode7 'clear Background operation
END SUB

'\$Page \$SubTitle:'SUB TestMode6'

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB TestMode6

SUB TestMode6

SHARED A2dBuff() AS INTEGER 'A/D data buffer
SHARED A2dChNum() AS INTEGER 'channel # buffer

DIM Op AS INTEGER, Status AS INTEGER, Count AS INTEGER, i AS
INTEGER

DasMode0 'init
DasMode1 0, 3 'scan 0 to 3
DasMode17 (100) 'sample @ 100 Hz
DasMode6 500, 0 'get 100 pts, only one cycle

DO

 DasMode8 Op, Status, Count
 PRINT Op, Status, Count
LOOP UNTIL Status = 0

DasMode9 500, 0 'Xfer to A2dBuff()

FOR i = 0 TO 499 STEP 4
 PRINT A2dChNum(i); A2dBuff(i), A2dChNum(i + 1); A2dBuff(i +
1), A2dChNum(i + 2); A2dBuff(i
+ 2), A2dChNum(i + 3); A2dBuff(i + 3)
NEXT i

DasMode7 'clear Background operation
END SUB

'\$Page \$SubTitle: 'SUB VentLoop'

SUBSTITUTE SHEET

160

2184992

VentLib.Bas - Ohio 7000 Ventilator Control

SUB VentLoop

SUB VentLoop

'Uses function to return avg of A/D channels when running
'in continuous acquisition mode 6.

'Modulates ventilator with data in common - if loaded.

SHARED VentArrNum AS INTEGER '# of items loaded in
VentArr()
SHARED VentDatFn\$ 'Data File loaded
SHARED VentArr() AS VentDat 'Buffer for Ventilator
modulation data

DIM Op AS INTEGER, Status AS INTEGER, Count AS INTEGER
DIM OldCount AS INTEGER, Done AS INTEGER, Chan AS INTEGER
DIM i AS INTEGER
DIM a\$, mask\$, mask1\$

DIM tStart AS SINGLE 'starting time
DIM tNow AS SINGLE 'current time

DIM NumCh AS INTEGER '# channels to acquire
DIM PtsPerCh AS INTEGER '# points per channel
DIM AvPer AS INTEGER 'Avg Period of A/D data

DIM VolMod, RateMod 'Vol & Rate Modulation level
DIM VolModV, RateModV 'Volts of Vol & Rate
Modulation level

DIM VolCtrl, RateCtrl 'Setting of ventilator Vol
& Rate Ctrl
DIM VolCtrlV, RateCtrlV 'Volts of ventilator Vol &
Rate Ctrl

SUBSTITUTE SHEET

2184992

161

VentLib.Bas - Ohio 7000 Ventilator Control

SUB VentLoop

DIM VolD2a, RateD2a
Rate

'D/A outputs for Volume and

DIM Vol, Rate
DIM VolV, RateV
Ventilator

'Current drive to Ventilator

'Current drive Voltage to

NumCh = 4
PtsPerCh = 20
AvPer = 10
DasMode0
DasMode7
operation

'channels 0 to 3

'acquire 20 pts per channel

'average previous 100 mS

'Init Driver

'Clear DMA or Interrupt

DasMode16 0, 0

'D/A's to zero

'Build Status Screen

COLOR 14, 1

CLS

a\$ = "Ventilation Control - Press ESC to Terminate"

LOCATE 1, 40 - LEN(a\$) / 2

PRINT a\$

IF VentArrNum > 0 THEN
if loaded

'Display Modulation file -

a\$ = "File: " + VentDatFn\$

LOCATE 25, 40 - LEN(a\$) / 2

PRINT a\$;

END IF

LOCATE 7, 30

PRINT "--Volume--"

LOCATE 8, 30

PRINT "L/Min Volt"

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB VentLoop

```
LOCATE 7, 45
PRINT "----Rate----"
LOCATE 8, 45
PRINT "Bth/M Volt"
LOCATE 10, 15
PRINT "BaseLine"
LOCATE 11, 15
PRINT "Modulation"
LOCATE 12, 15
PRINT "D/A Drive"
LOCATE 13, 15
PRINT "Output"
mask$ = "##.## #.###"
mask1$ = "      #.###"

DasModel 0, NumCh - 1           'Setup Das16 for channels 0-3
DasModel7 (100)                 '100 Hz/Ch sampling rate

'Use Mode6 to get 200 mS of Data - speed high - Background
DasMode6 NumCh * PtsPerCh, 1    'get 10 pts per ch, cycle
continuously

'Wait until filled A/D buffer once - 200 mS

DasMode8 Op, Status, OldCount   'look at A/D buffer index
pointer
Done = 0
DO
    DasMode8 Op, Status, Count   'look at A/D buffer index
pointer

IF Count < OldCount THEN        'wrapped around!
    Done = -1
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB VentLoop

ELSE

OldCount = Count

'wait until wrap around

END IF

LOOP UNTIL Done

VolMod = 2

'Zero Modulation Levels

RateMod = 6

'What to do

i = 0

DO

IF i = 0 THEN

'Init

i = 1

tStart = TIMER

'remember starting time

END IF

'If data loaded, see if it is time to change modulation level

IF VentArrNum > 0 THEN

tNow = TIMER - tStart

IF tNow < 0 THEN tNow = tNow + 86400 'working late eh!

IF tNow >= VentArr(i).Time THEN

'change modulation

' LOCATE 17, 30

' PRINT USING "####.# ! ##.##"; tNow; VentArr(i).Mode;

VentArr(i).Modulation

SELECT CASE VentArr(i).Mode

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB VentLoop

```
        CASE "R"
            RateMod = VentArr(i).Modulation
            i = i + 1
        CASE "V"
            VolMod = VentArr(i).Modulation
            i = i + 1
    END SELECT
END IF

'When we reach the end of the data send [VentArr()] - wrap
arround
    IF i > VentArrNum THEN i = 0

END IF

t1 = TIMER
IF 1 = 1 THEN                                'Volume

    'Read Volume control - avg of previous 100 mS
    DasMode8 Op, Status, Count                'determine current index
pointer
    DasMode9 NumCh * PtsPerCh, 0              'Xfer to A2dBuff(),
A2dChNum()

    'Read Volume control position
    Chan = 0
    VolCtrlV = A2dToVolt(AvgA2dChan$(Count, Chan, AvPer, NumCh,
PtsPerCh))
    VolCtrl = VolToV(1, VolCtrlV)

    'Read Volume Ventilator is running at
    Chan = 1
    VolV = A2dToVolt(AvgA2dChan$(Count, Chan, AvPer, NumCh,
PtsPerCh))
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB VentLoop

Vol = VolToV(1, VolV)

'Calculate Volume Modulation into volts

VolModV = VolToV(2, VolMod)

'Calculate D/A drive to Achieve Modulation Level

IF VolModV > VolCtrlV THEN 'Modulation > set point

VolD2a = VolModV - VolCtrlV

ELSE 'Modulation <= set point

VolD2a = 0

END IF

'Send the Volume Drive to the Ventilator via D/A channel

0

DasMode15 0, VoltToD2a*(VolD2a)

LOCATE 10, 30

PRINT USING mask\$; VolCtrl; VolCtrlV

LOCATE 11, 30

PRINT USING mask\$; VolMod; VolModV

LOCATE 12, 30

PRINT USING mask1\$; VolD2a

LOCATE 13, 30

PRINT USING mask\$; Vol; VolV

END IF

IF 2 = 2 THEN

'Rate

'Read Rate control - avg of previous 100 mS

DasMode8 Op, Status, Count

'determine current index

pointer

DasMode9 NumCh * PtsPerCh, 0

'Xfer to A2dBuff(),

A2dChNum()

SUBSTITUTE SHEET

2184992

166

VentLib.Bas - Ohio 7000 Ventilator Control
SUB VentLoop

```
'Read Rate control position
Chan = 2
RateCtrlV = A2dToVolt(AvgA2dChan$(Count, Chan, AvPer,
NumCh, PtsPerCh))
RateCtrl = RateToV(1, RateCtrlV)

'Read Rate Ventilator is running at
Chan = 3
RateV = A2dToVolt(AvgA2dChan$(Count, Chan, AvPer, NumCh,
PtsPerCh))
Rate = RateToV(1, RateV)

'Calculate Rate Modulation into volts
RateModV = RateToV(2, RateMod)

'Calculate D/A drive to Achieve Modulation Level
IF RateModV > RateCtrlV THEN 'Modulation > set point
    RateD2a = RateModV - RateCtrlV
ELSE 'Modulation <= set point
    RateD2a = 0
END IF

'Send the Rate Drive to the Ventilator via D/A channel 1
DasModel5 1, VoltToD2a$(RateD2a)

LOCATE 10, 45
PRINT USING mask$; RateCtrl; RateCtrlV

LOCATE 11, 45
PRINT USING mask$; RateMod; RateModV
LOCATE 12, 45
PRINT USING mask1$; RateD2a
```

SUBSTITUTE SHEET

2184992

167

VentLib.Bas - Ohio 7000 Ventilator Control
SUB VentLoop

LOCATE 13, 45
PRINT USING mask\$; Rate; RateV
END IF
t2 = TIMER

LOCATE 15, 30: PRINT USING "loop time = ##.##"; t2 - t1

a\$ = INKEY\$
LOOP UNTIL a\$ = CHR\$(27)

DasModel6 0, 0 'Zero D/A's
DasMode7 'Disable Interrupt

END SUB

'\$Page \$SubTitle:'SUB VentPlotGrid'

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB VentPlotGrid

SUB VentPlotGrid

 SHARED Vc AS Config, InitRows AS INTEGER, BestMode AS INTEGER,
Available AS STRING

 SHARED VentArrNum AS INTEGER '# of items loaded in
VentArr()

 SHARED VentDatFn\$ 'Data File loaded

 SHARED VentArr() AS VentDat 'Buffer for Ventilator
modulation data

 DIM txtX AS INTEGER, txtY AS INTEGER

 DIM x1 AS INTEGER, x2 AS INTEGER, y1 AS INTEGER, y2 AS INTEGER

 DIM i AS INTEGER

 DIM r AS SINGLE

 SCREEN BestMode

 SetConfig BestMode 'Fill in VC Parameters

 COLOR 1, 1 'White BackGround, Palette 1

 CLS

 'Calculate Char size in this mode

 txtX = (Vc.XPix + 1) \ Vc.TCol

 txtY = (Vc.YPix + 1) \ Vc.TRow

 'Label Header

 LOCATE 1, 7

 PRINT "File: "; VentDatFn\$;

 'Label Left Vertical Axis - Volume

 LOCATE 2, 2: PRINT "Vol";

 LOCATE 3, 3: PRINT "30";

 LOCATE Vc.TRow - 2, 3: PRINT " 0";

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control
SUB VentPlotGrid

```
'Tick Left Axis [0..30] every 5
x1 = 4 * txtX
x2 = 5 * txtX
y1 = 2 * txtY + (txtY) \ 2
y2 = (Vc.TRow - 2) * txtY - (txtY) \ 2
VIEW (x1, y1)-(x2, y2)
WINDOW (0, 1)-(1, 7)
FOR i = 1 TO 7
    LINE (0, i)-(1, i), 1
NEXT i

'Label Right Vertical Axis - Rate
LOCATE 2, Vc.TCol - 4: PRINT "Rate";
LOCATE 3, Vc.TCol - 4: PRINT "40";
LOCATE Vc.TRow - 2, Vc.TCol - 4: PRINT " 0";

'Tick Right Axis [0..40] every 5
x1 = (Vc.TCol - 5) * txtX
x2 = (Vc.TCol - 6) * txtX
y1 = 2 * txtY + (txtY) \ 2
y2 = (Vc.TRow - 2) * txtY - (txtY) \ 2
VIEW (x1, y1)-(x2, y2)
WINDOW (0, 1)-(1, 9)
FOR i = 1 TO 9
    LINE (0, i)-(1, i), 2
NEXT i

'Label X Axis
LOCATE Vc.TRow, (Vc.TCol \ 2) - 2
PRINT "(Sec)";
LOCATE Vc.TRow, 6
PRINT "0.0";
LOCATE Vc.TRow, Vc.TCol - 4 - 6
PRINT FUse$(VentArr(VentArrNum).Time, 6, 1);
```

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

SUB VentPlotGrid

'Tick Horizontal Axis by 100 Secs

x1 = 5 * txtX

x2 = (Vc.TCol - 6) * txtX

y1 = (Vc.TRow - 2) * txtY - (txtY) \ 2

y2 = y1 + (txtY \ 2)

VIEW (x1, y1)-(x2, y2)

WINDOW (0, 0)-(VentArr(VentArrNum).Time, 1)

FOR r = 0 TO VentArr(VentArrNum).Time STEP 100

LINE (r, 0)-(r, 1), 3

NEXT r

'Exit with the Plot Window defined

x1 = 5 * txtX' - 1

x2 = (Vc.TCol - 6) * txtX' + 1

y1 = 2 * txtY + (txtY) \ 2' - 1

y2 = (Vc.TRow - 2) * txtY - (txtY) \ 2' + 1

VIEW (x1, y1)-(x2, y2), 3, 3

END SUB

'\$Page \$SubTitle:'FUNCTION VolToV'

SUBSTITUTE SHEET

2184992

171

VentLib.Bas - Ohio 7000 Ventilator Control
FUNCTION VolToV

FUNCTION VolToV (Opt AS INTEGER, Volt)

'Sept 29/93

'-----
'Function Library for Ohio 7000 Ventilator Control
'-----

'Based on observations of Litre/Min Volume control

' Y = a + bX

' mV = 66.19015 + 66.71872 * Vol R^2=0.998923

' V = 6.61901478E-02 + 6.67187192E-02 * Vol

R^2=9.98922670E-01

' Vol Regression Output:

' Constant 6.61901478E-02

' Std Err of Y Est 1.89986311E-02

' R Squared 9.98922670E-01

' No. of Observations 29

' Degrees of Freedom 27

'

' X Coefficient(s) 6.67187192E-02

' Std Err of Coef. 4.21671541E-04

DIM V, Vol, a, b

' V = 6.61901478E-02 + 6.67187192E-02 * Vol

R^2=9.98922670E-01

a = .0661901478#

b = .0667187192#

SELECT CASE Opt%

CASE 1 'Volt to Vol

SUBSTITUTE SHEET

VentLib.Bas - Ohio 7000 Ventilator Control

FUNCTION VolToV

V = Volt

Vol = (V - a) / b ' x = (Y - a) / b

VolToV = Vol

CASE 2 'Vol to V

Vol = Volt

V = a + b * Vol ' y = a + bx

VolToV = V

END SELECT

END FUNCTION

43949 Bytes Available

28608 Bytes Free

0 Warning Error(s)

0 Severe Error(s)

SUBSTITUTE SHEET

Table III

Variable	Control Ventilator	Computer Ventilator	p-value
Weight (kg)	21.7 ± 2.8	23.4 ± 1.3	ns
Oleic Acid Infused (ml/kg)	0.20 ± 0.05	0.24 ± 0.11	ns
Mean Airway Pressure (cm H ₂ O)	12.02 ± 0.54	11.41 ± 0.39	ns
Mean Peak Airway Pressure (cm H ₂ O)	59.5 ± 1.3	56.6 ± 3.0	ns
Wet:Dry Weight Ratio	10.1 ± 1.1	9.2 ± 1.2	ns

Mean ± S.D.

Control Group n = 6
Computer Group n = 7
except for mean airway pressure and mean peak airway pressure
Control Group n = 4
Computer Group n = 3

Table IV

Variable	Baseline	Time (minutes) Following Oleic Acid Infusion						
		0	30	60	90	120	150	180
Temp (°C)	Computer	38.3±1.0*	38.4±1.1*	38.5±1.0	38.4±0.9*	38.5±1.0	38.6±1.1	38.6±1.3
	Control	38.0±1.2*	38.3±1.1*	38.4±1.2	38.5±1.0*	38.5±1.4	38.9±1.2	39.0±1.3
Temp (°C)	Computer	37.0±0.9	37.3±0.6	37.1±0.8	37.3±0.5	37.2±0.5	37.2±0.7	37.2±1.1
	Control	37.1±1.1*	37.5±0.9*	37.3±0.7	37.4±0.8*	37.7±0.7	37.7±1.0	37.9±1.0
Hgb (g%)	Computer	11.6±1*	11.8±.8*	11.9±.7*	11.6±.7*	11.7±.5*	11.8±.8*	12.3±.7*
	Control	12.6±1*+	13.0±1*+	13.3±1.3	13.7±1*+	13.9±1*+	14.1±1*+	14.5±.8*
pH	Computer	7.37±.04*	7.38±.05*	7.36±.07	7.38±.06*	7.37±.04	7.37±.04	7.34±.05
	Control	7.38±.03*	7.35±.02*	7.33±.02	7.33±.03*	7.33±.03	7.31±.06	7.29±.07

Mean ± S.D.
+ P < 0.05 between Groups
Bonferroni's correction applied for multiple comparisons

Computer Group n = 7
Control Group n = 6

Temp = blood temperature
Temp = nasopharyngeal temperature

2184992

Table V

Variable	Baseline	Time (minutes) Following Oleic Acid Infusion						
		0	30	60	90	120	150	180
MAP (mm Hg) Computer Control	91 ± 15	78 ±	75 ± 14*	75 ± 12*	73 ± 10*	75 ± 12*	76 ± 9*	77 ± 8*
	92 ± 14	11*+	69 ± 20*	71 ± 14*	71 ± 15*	78 ± 19*	79 ± 20*	78 ±
		68 ± 16*						17*
MPAP (mm Hg) Computer Control	22 ± 5	40 ± 4*	39 ± 3*	39 ± 5*	37 ± 5*	39 ± 6*	40 ± 7*	42 ± 4*
	19 ± 2+	37 ± 4*	37 ± 6*	40 ± 5*	40 ± 3*+	40 ± 5*	43 ± 2*	43 ± 3*
PCWP (mm Hg) Computer Control	10 ± 1	11 ± 1	11 ± .6	10 ± 1	10 ± .5	11 ± 1	12 ± 3	11 ± 3
	11 ± 1	12 ± 1	11 ± 2	11 ± 1	11 ± .8	10 ± 1	10 ± 2+	10 ± 1
PVR (mm Hg·l ⁻¹ ·min) Computer Control	3 ± 1	12 ± 1*	12 ± 2*	12 ± 2*	10 ± 2*	12 ± 4*	12 ± 5*	11 ± 4*
	2 ± .7	13 ± .4*	13 ± 6*	15 ± 6*	14 ± 5*+	15 ± 5*+	15 ± 7*+	15 ±
								5*+
CO (l·min ⁻¹) Computer Control	4.2 ± .3	2.5 ±	2.6 ±	2.5 ±	2.6 ± .3*	2.5 ±	2.5 ±	2.7 ±
	4.0 ± .9	.4*	.4*	.4*	2.3 ± .6*	.3*	.4*	.4*
		2.4 ±	2.5 ±	2.5 ±		2.5 ±	2.5 ±	2.5 ±
		.6*	.9*	.9*		.8*	1.0*	.8*

Mean ± S.D.

* P < 0.05 within Groups

+ P < 0.05 between Groups

Bonferroni's correction applied for multiple comparisons

Computer Group n = 7

Control Group n = 6

MPAP - Mean Pulmonary Artery Pressure

PCWP - Pulmonary Capillary Wedge Pressure

PVR - Pulmonary Vascular Resistance

CO - Cardiac Output

SUBSTITUTE SHEET

Table VI

Variable	Baseline	Time (minutes) Following Oleic Acid Infusion						
		0	30	60	90	120	150	180
PaCO ₂ (mm Hg) Computer Control	17.5 ±	13.7 ±	14.4 ±	16.1 ±	14.7 ±	15.0 ±	13.1 ±	15.2 ± 2.7
	4.8	2.6	2.5	3.2	2.9	2.2	2.3	17.3 ± 6.1
	19.8 ±	16.1 ±	14.5 ±	15.4 ±	15.7 ±	16.8 ±	17.8 ±	
	3.7	4.1	2.5	3.6	3.8	4.9	6.8	
PaCO ₂ (mm Hg) Computer Control	36.3 ± 2	46.6 ± 3*	45.3 ± 6*	47.7 ± 8*	45.8 ± 8*	46.0 ± 4*	46.8 ± 6*	49.9 ± 6*
	36.9 ± 2	43.6 ± 1*	46.1 ± 4*	48.2 ± 7*	48.1 ± 3*	49.0 ± 5*	49.0 ± 7*	51.2 ± 13*
PaO ₂ (mm Hg) Computer Control	558 ± 31	133 ± 39*	140 ± 68*	158 ± 93*	203 ± 88*	197 ±	162 ±	116 ± 64*
	556 ± 57	112 ± 56*	103 ± 56*	92 ± 51**	90 ± 53**	124*	119*	65 ± 15*
QS/QT Computer Control	10.0 ±	18.4 ±	17.4 ±	16.8 ±	15.9 ±	15.9 ±	16.7 ±	17.8 ±
	4.6	3.9*	4.0*	4.1*	3.5*	4.2*	3.4*	3.9*
	9.0 ± 3.0	18.6 ±	19.0 ±	17.7 ±	17.9 ±	18.3 ±	17.6 ±	18.3 ±
		4.6*	5.7*	3.3*	4.1*	4.0*	4.6*	5.1*
VD/VT Computer Control	52.1 ±	70.5 ±	68.0 ±	66.0 ±	67.5 ±	67.3 ±	71.6 ±	69.6 ±
	10.2	5.8*	5.2	5.0*	6.1*	5.2*	5.3*	4.1*
	45.6 ±	63.3 ±	68.3 ±	67.9 ±	67.5 ±	65.6 ±	64.4 ±	66.9 ±
	12.1	8.5*	6.0*	7.4*	7.2*	8.1*	10*	4.6*

Mean ± S.D.

* P < 0.05 within Groups

+ P < 0.05 between Groups

Bonferroni's correction applied for multiple comparisons

Computer Group n = 7

Control Group n = 6

PaCO₂ = end expired CO₂

QS/QT = shunt fraction

VD/VT = dead space ventilation

SUBSTITUTE SHEET

Table VII

Variable	Hypothermia	Rewarming
Temp °C		
Computer	28.4 ± .3	35.6 ± .82*
Control	28.0 ± .24	35.0 ± .91*
MAP (mm Hg)		
Computer	91 ± 26	93 ± 17
Control	81 ± 13	82 ± 10+
CSFP (mm Hg)		
Computer	9.1 ± 3	14.6 ± 3.7*
Control	8.2 ± 4.3	12.3 ± 4.1*+
CPP (mm Hg)		
Computer	82 ± 26	80 ± 19
Control	72 ± 15	69 ± 11

Computer n=6

Control n=6

* p < 0.05 within groups

+ p < 0.05 between groups

SUBSTITUTE SHEET

178 2184992

Table VIII

Variable	Hypothermia	Rewarming
Hgb (g/dl)		
Computer	7.4 ± .7	8.1 ± .9
Control	7.7 ± 1.1	7.7 ± 1
PaCO ₂ (mm Hg)		
Computer	37 ± 2	36 ± 1
Control	38 ± 2	36 ± 6
pH		
Computer	7.35 ± .01	7.35 ± .02
Control	7.33 ± .02	7.35 ± .08
Cont Diff (Vol %)		
Computer	4.4 ± 1.2	4.3 ± 1.1
Control	3.6 ± 0.7	4.9 ± 1.1*
SSS PO ₂ (mm Hg)		
Computer	44 ± 6	42 ± 7*
Control	44 ± 3	38 ± 4*+
SSS Sat (%)		
Computer	68 ± 11	67 ± 10
Control	69 ± 6	60 ± 7*+

Computer n = 6

Control n = 6

* p < 0.05 within groups

+ p < 0.05 between groups

SUBSTITUTE SHEET

Table IX

Variable	Hypothermia	Rewarming
tCBF (ml.g ⁻¹ .min ⁻¹)		
Computer	.18 ± .08	.36 ± .07*
Control	.17 ± .06	.33 ± .06*
hCBF (ml.g ⁻¹ .min ⁻¹)		
Computer	.17 ± .1	.35 ± .08*
Control	.17 ± .06	.32 ± .06*
bsCBF (ml.g ⁻¹ .min ⁻¹)		
Computer	.21 ± .11	.41 ± .07*
Control	.20 ± .07	.38 ± .08*
FLOW:METABOLISM		
Computer	23.1 ± 5.2	23.9 ± 5.0
Control	28.6 ± 5.6	22.3 ± 4.3*
CMRO ₂		
Computer	.008 ± .004	.016 ± .008*
Control	.006 ± .002	.014 ± .002*

Computer n = 6

Control n = 6

* p < 0.005 within groups

+ p < 0.005 between groups

SUBSTITUTE SHEET

The embodiments of the invention in which an exclusive property or privilege is claimed are defined as follows:

1. Apparatus for controlling the flow of a biological fluid to an organ, which comprises:

means for establishing a predetermined pattern of variations over time of instantaneous changes in flow of a biological fluid to an independently-functioning normal organ of a mammalian species,

means for generating a variable control parameter for regulation of flow of the biological fluid to an organ during controlled life support conditions in accordance with the predetermined pattern, and

means for controlling the flow of the biological fluid to the organ during controlled life support conditions in accordance with the variable control parameter.

2. The apparatus of claim 1 for use in controlling the flow of blood by a pump to a body during cardiopulmonary by-pass wherein said means for establishing a predetermined pattern comprises means for establishing a predetermined pattern of variation over time of instantaneous blood pressure and heart rate of an independently-functioning heart of a mammalian species; said means for generating a variable control parameter comprises means for generating a signal corresponding in value to an individually-determined blood pressure for a period of time corresponding to the heart rate for a difference between such individually-determined blood pressure and the next individually-determined blood pressure of said predetermined pattern and means for

generating a control voltage corresponding in magnitude to said signal; and said means for controlling the flow of the biological fluid comprises means for applying said control voltage to said pump to provide an output of blood from said pump to the body during cardiopulmonary bypass of a pressure proportional to the magnitude of the signal for said period of time and means for repeating the steps of generating a signal, generating a control voltage and applying the control voltage to the pump for each next individually-determined blood pressure of said predetermined pattern.

3. The apparatus of claim 1 for use in controlling the flow of ventilating gas from a ventilator device to the lungs of a body during controlled life support conditions wherein said means for establishing a predetermined pattern comprises means for establishing a predetermined pattern of variation over time of instantaneous respiratory rate and tidal volume of the independently-functioning healthy lungs of a mammalian species; said means for generating a variable control parameter comprises means for generating a signal corresponding in value to an individually-determined respiratory rate and tidal volume in said predetermined pattern and means for generating control voltage corresponding in magnitude to said signal; and said means for controlling the flow of biological fluid comprises means for applying the control voltage to the ventilator device to provide an output of ventilating gas from the ventilating device of a respiratory rate proportional to the said signal and means for repeating the steps of generating a signal, generating a control voltage and

applying the control voltage to the ventilator device for each next individually-determined respiratory rate of the predetermined pattern.




FIG.1.

SUBSTITUTE SHEET

2184992

2/39

SYSTOLIC BLOOD PRESSURE VARIABILITY
Input File for Computer-controller




FIG. 2.

SUBSTITUTE SHEET

3/39

2184992




FIG. 3.

SUBSTITUTE SHEET

4/39

2184992




FIG.4.

SUBSTITUTE SHEET

5/39

2184992




FIG. 5.

SUBSTITUTE SHEET

6/39

2184992




FIG. 6.

SUBSTITUTE SHEET

7/39

2184992




FIG. 7.

SUBSTITUTE SHEET

8/39

2184992





FIG.8.

SUBSTITUTE SHEET

9/39

2184992



Seconds

FIG. 9.

SUBSTITUTE SHEET

10/39

2184992




FIG.10.

SUBSTITUTE SHEET

2184992

11/39





FIG.11.

SUBSTITUTE SHEET

2184992

12/39



SUBSTITUTE SHEET

2184992

13/39





FIG.12B.

SUBSTITUTE SHEET

2184992

14/39



SUBSTITUTE SHEET

FIG.12C.

2184992

15/39





FIG.12D.

16/39

2184992



SUBSTITUTE SHEET

FIG.12E.

2184992

17/39




FIG.13 A.

SUBSTITUTE SHEET

18/39






FIG.13B.

SUBSTITUTE SHEET

2184992





FIG.13C.

SUBSTITUTE SHEET

20/39

2184992



- OUTPUT RANGE IS A MAX OF 2 VOLTS

- REGISTOR VALUES ON CONTROLS
MAY BE INCORRECT

- ONLY SIGNAL ROUTING OF
RATE, VOL, & IE CORRECT

FIG.13D.

SUBSTITUTE SHEET

2184992

21/39




FIG.13E.

SUBSTITUTE SHEET




FIG. 13F.

23/39 2184992




FIG.14 A.

SUBSTITUTE SHEET




FIG.14B.

SUBSTITUTE SHEET

2184992

25/39





FIG.14C.

SUBSTITUTE SHEET

2184992

26/39



VOL CTRL

FIG.14D.

SUBSTITUTE SHEET

27/39

2184992




FIG.14E.

SUBSTITUTE SHEET

28/39

2184992




FIG.14F.

SUBSTITUTE SHEET




FIG.15A.

SUBSTITUTE SHEET




FIG.15B.

SUBSTITUTE SHEET




FIG.15C.

2184992

32/39




FIG.15D.

SUBSTITUTE SHEET

33/39




FIG.15E.

SUBSTITUTE SHEET

34/39





FIG.15F.

SUBSTITUTE SHEET

2184992

35/39



DAS16 'JUMPER BOX'

FIG.16.

SUBSTITUTE SHEET

2184992

36/39




FIG.17 A.

SUBSTITUTE SHEET

2184992

37/39




FIG.17B.

SUBSTITUTE SHEET

2184992

38/39





FIG.17C.

TIE CHAN 8-15 TO GND

SUBSTITUTE SHEET

FIG.17D.



SUBSTITUTE SHEET

OBSERVE &
RECORD
SYSTOLIC
PRESSURE

LOAD INTO
COMPUTER
MEMORY

CALCULATE
MIN, MAX AVG
SYSTOLIC PRESSURE

PROCESS PEAK
PRESSURES

SET MIN &
MAX RATES

D/A
CONVERTOR

CONTROL
VOLTAGE

PUMP

COMPUTER
OPERATIONS

