



(51) International Patent Classification:

A61K 9/08 (2006.01) A61P 7/04 (2006.01)

A61K 48/00 (2006.01) A61K 47/02 (2006.01)

A61K 47/26 (2006.01) A61K 47/10 (2017.01)

(21) International Application Number:

PCT/US2021/064222

(22) International Filing Date:

17 December 2021 (17.12.2021)

(25) Filing Language:

English

(26) Publication Language:

English

(30) Priority Data:

63/127,826 18 December 2020 (18.12.2020) US

(71) Applicants: SANGAMO THERAPEUTICS, INC.

[US/US]; 7000 Marina Boulevard, Brisbane, California 94005 (US). **PFIZER INC.** [US/US]; 235 East 42nd Street, New York, New York 10017 (US).

(72) Inventors: CONNER, Jessica Eileen; c/o Pfizer Inc.,

875 Chesterfield Parkway West, Chesterfield, Missouri 63017 (US). **CRAWFORD, Lindsey Anne**; c/o Pfizer Inc., 1 Burtt Road, Andover, Massachusetts 01801 (US). **DAMITZ, Robert**; c/o Sangamo Therapeutics, Inc., 7000 Marina Boulevard, Brisbane, California 94005 (US). **DAVIS, Brendan Michael**; c/o Pfizer Inc., 1 Burtt Road, Andover, Massachusetts 01810 (US). **HODGE, Cody Michael**; c/o Pfizer Inc., 875 Chesterfield Parkway West, Chesterfield, Missouri 63017 (US). **KIMMEL II, Michael Leland**; c/o Pfizer Inc., 875 Chesterfield Parkway West, Chesterfield, Missouri 63017 (US). **QURESHI WIL-LARD, Tihami**; c/o Pfizer Inc., 875 Chesterfield Parkway West, Chesterfield, Missouri 63017 (US). **RAMSEY, Phillip**; c/o Sangamo Therapeutics, Inc., 7000 Marina Boulevard, Brisbane, California 94005 (US). **THORNE, Daniel Joseph**; c/o Pfizer Inc., 875 Chesterfield Parkway West, Chesterfield, Missouri 63017 (US). **YOUNG, Anthony Lee**; c/o Pfizer Inc., 875 Chesterfield Parkway West, Chesterfield, Missouri 63017 (US).

(74) Agent: LI, Z. Ying et al.; Steptoe & Johnson LLP, 1330

Connecticut Avenue, N.W., Washington, District of Columbia 20036 (US).

(81) Designated States (unless otherwise indicated, for every

kind of national protection available): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DJ, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IR, IS, IT, JO, JP, KE, KG, KH, KN, KP, KR, KW, KZ, LA, LC, LK, LR, LS, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO,

NZ, OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SA, SC, SD, SE, SG, SK, SL, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, WS, ZA, ZM, ZW.

(84) Designated States (unless otherwise indicated, for every

kind of regional protection available): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, ST, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, KM, ML, MR, NE, SN, TD, TG).

Declarations under Rule 4.17:

- as to applicant's entitlement to apply for and be granted a patent (Rule 4.17(ii))
- as to the applicant's entitlement to claim the priority of the earlier application (Rule 4.17(iii))
- of inventorship (Rule 4.17(iv))

Published:

- with international search report (Art. 21(3))
- with sequence listing part of description (Rule 5.2(a))

(54) Title: IMPROVED PHARMACEUTICAL COMPOSITIONS CONTAINING ADENO-ASSOCIATED VIRAL VECTOR

(57) Abstract: The present invention provides compositions comprising a recombinant AAV and one or more pharmaceutically acceptable excipients. The compositions have improved stability and shelf life as compared to other AAV compositions.

IMPROVED PHARMACEUTICAL COMPOSITIONS CONTAINING ADENO-ASSOCIATED VIRAL VECTOR

CROSS-REFERENCE TO RELATED APPLICATION

[0001] This application claims priority from U.S. Application No. 63/127,826, filed December 18, 2020, the content of which is incorporated by reference herein in its entirety.

SEQUENCE LISTING

[0002] The instant application contains a Sequence Listing which has been submitted electronically in ASCII format and is hereby incorporated by reference in its entirety. Said ASCII copy, created on December 1, 2021, is named 025297_WO029_SL.txt and is 19,967 bytes in size.

BACKGROUND OF THE INVENTION

[0003] Gene therapy is a promising approach to treat genetic diseases. It introduces a healthy copy of the defective gene to the patient, or inactivates a mutated gene that is functioning improperly. A genetic disease that is of particular interest for gene therapy research is hemophilia A. Hemophilia A, also called classic hemophilia, is a X-linked genetic disease in which the blood clotting process is impaired due to a missing or defective gene encoding factor VIII. Patients with hemophilia A can bleed internally (e.g., into joints and muscles) or externally (e.g., due to minor cuts, trauma, or dental procedures). Normal plasma levels of factor VIII range from 50% to 100%. Mild hemophilia A is characterized by a level of 6% to 49% of factor VIII in the blood; patients typically bleed only after serious injury, trauma, or surgery. Moderate hemophilia A is characterized by a level of 1% to 5% of factor VIII in the blood; patients have bleeding episodes after injuries. Severe hemophilia A is characterized by a level of less than 1% of factor VIII in the blood; patients experience bleeding after injuries and also have frequent spontaneous bleeding episodes, often in their joints and muscles. *See also* the website of the National Hemophilia Foundation. Hemophilia A is commonly treated on-demand (upon

bleeding) or prophylactically with replacement factor VIII. Some patients develop alloantibodies (aka inhibitors) against the replacement factor, rendering the therapy ineffective. Current therapies are burdensome because they require frequent intravenous injections and also are limited by resources in developing countries. Thus, gene therapy offers a promising approach to treating hemophilia A.

[0004] Recombinant adeno-associated virus (rAAV) has been explored as a platform for gene delivery in gene therapy. Adeno-associated virus (AAV) is a small non-enveloped virus belonging to the family *Parvoviridae* and the genus *Dependoparvovirus*. The virus is composed of a single-stranded DNA genome packaged into capsids assembled from three capsid proteins – viral protein (VP) 1, VP2, and VP3. Formulating rAAV preparations into pharmaceutical compositions for clinical use has been a challenge. Commonly used rAAV formulations have been observed to form visible precipitates over time or during lab-simulated stress tests. Left unaddressed, these precipitates represent a risk to patient safety. Physical stability issues observed with some rAAV formulations also could negatively impact the efficacy of products over the shelf life required for storage, shipment, and administration to patients (*see, e.g., Wright et al., Molecular Therapy* (2005) 12(1):171-8; Croyle et al., *Gene Therapy* (2001) 8:1281-90). Thus, there is a need for developing improved formulations for rAAV vectors so that their therapeutic potential in gene therapy can be fully implemented.

SUMMARY OF THE INVENTION

[0005] The present disclosure provides stable rAAV vector formulations suitable for clinical administration. In an aspect, the disclosure provides a pharmaceutical composition comprising an rAAV vector, sodium chloride (NaCl), potassium chloride (KCl), disodium phosphate (Na_2HPO_4), monopotassium phosphate (KH_2PO_4), magnesium chloride (MgCl_2), a polyol (e.g., sucrose), and a poloxamer (e.g., poloxamer 188), optionally wherein the composition comprises no more than about 0.1 mM calcium chloride and has a pH of about 7.1 to about 7.5.

[0006] In some embodiments, the composition contains about 0.1 to about 2.0 mM (e.g., about 0.5 mM or higher, about 1.3 mM or higher, or about 1.4 mM) magnesium chloride.

[0007] In some embodiments, the composition contains about 150 to about 200 mM, optionally about 172 mM, sodium chloride.

[0008] In some embodiments, the composition contains about 2.5 to about 3.0 mM, optionally about 2.7 mM, potassium chloride.

[0009] In some embodiments, the composition contains about 5 to about 10 mM, optionally about 8 mM, disodium phosphate.

[0010] In some embodiments, the composition contains about 1.0 to about 2.0 mM, optionally about 1.5 mM, monopotassium phosphate.

[0011] In some embodiments, the composition contains about 0.5% to about 2% (w/v), optionally about 1% (w/v), sucrose.

[0012] In some embodiments, the composition contains about 0.01% to about 0.1% (w/v), optionally about 0.05% (w/v), poloxamer 188.

[0013] In particular embodiments, the present disclosure provides a pharmaceutical composition comprising an rAAV vector, about 171.81 mM sodium chloride, about 2.68 mM potassium chloride, about 8.10 mM disodium phosphate, about 1.47 mM monopotassium phosphate, about 1.40 mM magnesium chloride, about 1.00% (w/v) sucrose, and about 0.05% (w/v) poloxamer 188, optionally wherein the composition comprises no more than about 0.1 mM calcium chloride and has a pH of about 7.1 to about 7.5.

[0014] In other particular embodiments, the present disclosure provides a pharmaceutical composition comprising an rAAV vector, about 172 mM sodium chloride, about 2.68 mM potassium chloride, about 8.10 mM disodium phosphate, about 1.47 mM monopotassium phosphate, about 0.49 mM magnesium chloride, about 1.00% (w/v) sucrose, and about 0.05% (w/v) poloxamer 188, optionally wherein the composition comprises no more than about 0.1 mM calcium chloride and has a pH of about 7.1 to about 7.5.

[0015] In some embodiments, the rAAV in the present compositions comprises a genome comprising an expression cassette for a therapeutic protein, such as a human factor VIII polypeptide (e.g., SEQ ID NO:1). In particular embodiments, the AAV genome comprises SEQ ID NO:2 or nucleotides 131-5,024 of SEQ ID NO:2.

[0016] In some embodiments, the composition contains the rAAV at about 1.0E+12 to about 1.0E+14 vector genomes (vg) per mL, optionally about 1.0E+13 to about 5.0E+13 vg per mL (e.g., about 1.0E+13 vg per mL).

[0017] In some embodiments, the rAAV comprises an AAV6 capsid protein (e.g., having an AAV6 capsid). In some embodiments, the AAV genome comprises inverted terminal repeats (ITR) from AAV2.

[0018] In another aspect, the present disclosure provides a vial comprising 5-10 mL, optionally 6.4 mL, of the present composition. The vial may be made of, e.g., cyclo-olefin copolymer, and/or may have an in-place thermoplastic elastomer stopper.

[0019] In another aspect, the present disclosure provides a method of treating a patient in need of a therapeutic protein, comprising administering to the patient the present composition. In some embodiments, the present disclosure provides a method of increasing the serum level of factor VIII in a human subject in need thereof (e.g., a human subject having hemophilia A), comprising administering intravenously to the human subject the present composition where the rAAV encodes a human factor VIII polypeptide. Also provided are the pharmaceutical compositions for use in such treatment methods, and the use of the compositions for the manufacture of a medicament for use in such methods.

[0020] Other features, objects, and advantages of the invention are apparent in the detailed description that follows. It should be understood, however, that the detailed description, while indicating embodiments and aspects of the invention, is given by way of illustration only, not limitation. Various changes and modification within the scope of the invention will become apparent to those skilled in the art from the detailed description.

BRIEF DESCRIPTION OF THE FIGURES

[0021] **FIGs. 1A-F** summarize the characterization of the particulate matter found in samples of formulations containing SB-525.

[0022] **FIG. 2** shows a schematic of buffer exchange by tangential flow filtration (TFF).

[0023] **FIG. 3** shows the formulation, fill, and finish steps for the AAV6 vector formulations.

[0024] **FIG. 4A** shows the vector genome (vg) titer and infectious titer of AAV6 vector formulations following multiple freeze thaw cycles. Error bars displayed are either $\pm$ the percentage RSD from $N \geq 3$ sample measurements for particle concentration by MADLS, or $\pm$ the accepted test variability for vg titer, median tissue culture infectious dose (TCID₅₀), ELISA, and SE-LC. RSD: relative standard deviation. MADLS: multi-angle dynamic light scattering. SE-LC: size exclusion liquid chromatography.

[0025] FIG. 4B shows the capsid titer and particle concentration of AAV6 vector formulations following multiple freeze/thaw cycles. Error bars displayed are either $\pm$ the percentage RSD from $N \geq 3$ sample measurements for particle concentration by MADLS, or $\pm$ the accepted test variability for vg titer, TCID₅₀, ELISA, and SE-LC.

[0026] FIG. 4C is a table showing the general and purity quality attribute results of freeze/thaw AAV6 vector formulation samples.

[0027] FIG. 4D is a table showing the strength quality attribute results of freeze/thaw AAV6 vector formulation samples.

[0028] FIG. 5A shows vg and infectious titer over up to 6 months at ambient (25°C/60% RH) conditions. RH: relative humidity.

[0029] FIG. 5B shows capsid titer and particle concentration over up to 6 months at ambient (25°C/60% RH) conditions.

[0030] FIG. 5C is a table showing the general and purity quality attribute results of ambient (25°C/60% RH) incubated samples. ND: no data. NAA: no analysis available.

[0031] FIG. 5D is a table showing the strength quality attribute results of ambient (25°C/60% RH) incubated samples.

[0032] FIG. 6A shows vg and infectious titer over 3 months at stressed (40°C/75% RH) conditions.

[0033] FIG. 6B shows capsid titer and particle concentration over up to 7 months at stressed (40°C/75% RH) conditions.

[0034] FIG. 6C is a table showing the general and purity quality attribute results of stressed (40°C/75%RH) incubated samples.

[0035] FIG. 6D is a table showing the strength quality attribute results of stressed (40°C/75% RH) incubated samples.

[0036] FIG. 7 is a table showing quality attribute with associated and failure criteria for stability samples.

[0037] FIG. 8 is a table showing the inter-sample formulation assessment summary.

[0038] FIG. 9 is a table showing the study endpoint trendline and formulation assessments for general and purity quality attributes.

[0039] FIG. 10 is a table showing the long term (24 month) stability of the SB-525 drug product at the intended storage temperature (-70°C).

DETAILED DESCRIPTION OF THE INVENTION

[0040] The present disclosure provides pharmaceutical compositions comprising AAV vectors and one or more pharmaceutically acceptable excipients. The present AAV vector compositions may comprise rAAV whose genome carries an expression cassette for a protein of interest (e.g., a therapeutic protein). The present inventors have unexpectedly discovered that AAV vector formulations that are essentially calcium free (e.g., having no added calcium in the formulation) have improved stability and shelf life as compared to prior composition. Calcium is typically included in prior formulations due to the conventional belief that calcium would improve the stability of AAV compositions (*see, e.g.,* Turnbull et al., *Hum Gene Ther.* (2000) 11(4):629-35; Cotmore et al., *J Virol.* (2010) 84(4):1945-56). The inventors have discovered that the present AAV vector compositions have improved appearance (e.g., clarity and lack of color), more stable pH, and less aggregates (as determined by product quality attributes under freeze/thaw cycles and accelerated stability conditions). Because the AAV vector product formulated without calcium has demonstrated better product stability, calcium ions are not required for product performance or stability.

I. Preparation of Recombinant AAV

[0041] The viral preparations described herein may be obtained by any known production systems, such as mammalian cell AAV production systems (e.g., those based on 293T or HEK293 cells) and insect cell AAV production systems (e.g., those based on sf9 insect cells and/or those using baculoviral helper vectors). The viral preparations may be purified from the cell cultures by using well known techniques such as discontinuous cesium chloride density gradients (*see, e.g.,* Grieger, *Mol Ther Methods Clin Dev.* (2016) 3:16002).

[0042] The present compositions may comprise AAV of any or a combination of a variety of AAV serotypes, such as AAV1, AAV2, AAV3, AAV3B, AAV4, AAV5, AAV6, AAV7, AAV8, AAV8.2, AAV9, AAVrh10, AAV10, and AAV11, as well as variants, hybrids, chimera or pseudo-types thereof. By “pseudo-typed” or “cross-packaged” rAAV is meant a recombinant AAV whose capsid is replaced with the capsid of another AAV serotype, to, for example, alter transduction efficacy or tropism profiles of the virus (*see, e.g.,* Balaji et al., *J Surg Res.* (2013) 184(1):691-8). By “chimeric” or “hybrid” rAAV is meant a recombinant AAV whose capsid is

assembled from capsid proteins derived from different serotypes and/or whose capsid proteins are chimeric proteins with sequences derived from different serotypes (e.g., serotypes 1 and 2; *see, e.g.*, Hauck et al., *Mol Ther.* (2003) 7(3):419-25). For example, the present compositions may comprise recombinant AAV whose genome such as the ITRs is derived from one serotype such as AAV2 while the capsids are derived from another serotype; e.g., AAV2/8, AAV2/5, AAV2/6, AAV2/9, or AAV2/6/9. *See, e.g.*, U.S. Pats. 7,198,951 and 9,585,971.

II. Formulation of Recombinant AAV

[0043] Once purified, the AAV preparations can be formulated as described herein, for example, by buffer exchange through tangential flow filtration, normal flow filtration using stir-cells, gel filtration, dialysis, column chromatography, and/or desalting columns, to arrive at a composition comprising the desired ingredients. By way of illustration, the purified viral preparation may be concentrated first by ultrafiltration (UF) and then diafiltrated (DF) with 10 times or more equivalent volumes of the desired aqueous formulation solution. *See also* the Examples below.

[0044] The formulation solution may comprise tonicity agents, stabilizing agents, surfactants, and buffering agents. Buffering agents may include, for example, acetate, succinate (e.g., disodium succinate hexahydrate), succinic acid, gluconate, citrate, histidine, acetic acid, phosphate, phosphoric acid, ascorbate, ascorbic acid, tartaric acid, malate, maleic acid, glycine, lactate, lactic acid, bicarbonate, carbonic acid, sodium benzoate, benzoic acid, edetate, imidazole, tris, and mixtures thereof. In some embodiments, the formulation solution contains sodium chloride and/or potassium chloride at, e.g., about 150-200 mM (e.g., about 150 mM, about 155 mM, about 160 mM, about 165 mM, about 168 mM, about 170 mM, about 171 mM, about 171.1 mM, about 171.2 mM, about 171.3 mM, about 171.4 mM, about 171.5 mM, about 171.6 mM, about 171.7 mM, about 171.8 mM, about 171.9 mM, or about 172 mM) and about 2.5-3.0 mM (e.g., about 2.5 mM, about 2.6 mM, about 2.61 mM, about 2.61 mM, about 2.63 mM, about 2.64 mM, about 2.65 mM, about 2.66 mM, about 2.67 mM, about 2.68 mM, about 2.69 mM, or about 2.7 mM, about 2.8 mM, about 2.9 mM, about 3.0 mM), respectively.

[0045] The formulation solution may be phosphate-buffered, e.g., by disodium phosphate and/or monopotassium phosphate. In some embodiments, the total phosphate ion concentration in the formulation solution is about 8-12 mM (e.g., about 9.6 mM or about 9.57 mM). In certain

embodiments, the formulation solution contains about 5-10 mM (e.g., about 5 mM, about 6 mM, about 7 mM, about 7.9 mM, about 8 mM or about 8.1 mM, about 8.2 mM, about 8.5 mM, about 9 mM, or about 10 mM) disodium phosphate and about 1-2 mM (e.g., about 1 mM, about 1.2 mM, about 1.3 mM, about 1.45 mM, about 1.47 mM, about 1.48 mM, about 1.49 mM, about 1.5, about 1.6 mM, about 1.7 mM, about 1.8 mM, about 1.9 Mm, or about 2 mM) monopotassium phosphate. The formulation may have a pH of about 6.5-8.0 (e.g., about 7.1-7.5, about 7.1, about 7.2, about 7.3, about 7.4, or about 7.5).

[0046] The formulation solution may contain magnesium but no added calcium. In some embodiments, the formulation solution contains about 0.1 to 2.0 mM (e.g., about 0.5 to about 1.4 mM, about 0.1 mM, about 0.2 mM, about 0.3 mM, about 0.4 mM, about 0.42 mM, about 0.44 mM, about 0.45 mM, about 0.46 mM, about 0.47 mM, about 0.48 mM, about 0.49 mM, about 0.5, about 0.55 mM, or about 0.6 mM,) magnesium chloride. Although the formulation solution contains no added calcium, the pharmaceutical composition made up from the AAV preparation and the formulation solution may have a trace amount of calcium that is carried over for the AAV manufacturing and purification process. For example, the pharmaceutical composition may contain no more than about 0.10 mM (e.g., no more than about 0.09, 0.07, 0.05, 0.03, or 0.01 mM) calcium when measured by a colorimetric assay. In some embodiments, the pharmaceutical composition contains no detectable calcium as measured by a colorimetric assay. In some embodiments, the pharmaceutical composition contains no calcium (i.e., 0 mM calcium).

[0047] The formulation solution may contain a polyol such as mannitol, trehalose, sorbitol, erythritol, isomalt, lactitol, maltitol, xylitol, glycerol, lactitol, ethylene glycol, propylene glycol, polyethylene glycol, inositol, fructose, glucose, mannose, sucrose, sorbose, xylose, lactose, maltose, dextran, pullulan, dextrin, cyclodextrins, soluble starch, hydroxyethyl starch, water-soluble glucans, or mixtures thereof. In some embodiments, the formulation solution contains about 0.5% to 2% (e.g., about 0.5%, about 0.6%, about 0.7%, about 0.8%, about 0.9%, or about 1%) (w/v) sucrose.

[0048] The formulation solution may contain a nonionic or ionic hydrophilic surfactant. Example of surfactants are a polysorbate, poloxamer, triton, sodium dodecyl sulfate, sodium laurel sulfate, sodium octyl glycoside, lauryl-sulfobetaine, myristyl-sulfobetaine, linoleyl-sulfobetaine, stearyl-sulfobetaine, lauryl-sarcosine, myristyl-sarcosine, linoleyl-sarcosine,

stearyl-sarcosine, linoleyl-betaine, myristyl-betaine, cetyl-betaine, lauroamidopropyl-betaine, cocamidopropyl-betaine, linoleamidopropyl-betaine, myristamidopropyl-betaine, palmidopropyl-betaine, isostearamidopropyl-betaine, myristamidopropyl-dimethylamine, palmidopropyl-dimethylamine, isostearamidopropyl-dimethylamine, sodium methyl cocoyl-aurate, disodium methyl oleyl-aurate, dihydroxypropyl PEG-5 linoleammonium chloride, polyethylene glycol, polypropylene glycol, sorbitan monoesters (e.g., Spans), acid esters of glycerol, and mixtures thereof. In some embodiments, the surfactant can be polysorbate (PS) 20, PS-21, PS-40, PS-60, PS-61, PS-65, PS-80, PS-81, PS-85, PEG-3350, poloxamer 188, and mixtures thereof. In some embodiments, the formulation solution contains about 0.01% to 0.1% (e.g., about 0.01%, about 0.02%, about 0.03%, about 0.04%, about 0.05%, about 0.06%, about 0.07%, about 0.08%, about 0.09%, or about 0.1%) (w/v) poloxamer 188.

[0049] In particular embodiments, the pharmaceutical composition comprises an AAV vector in Formulation F2 or F3. The ingredients for F2 and F3 are shown in **Table A** below.

Dulbecco's phosphate buffered saline (DPBS) with calcium and magnesium is a commonly used formulation for cell culture applications. F0 is another prior, calcium-containing formulation.

As shown in the Examples below, F2 and F3 are superior to DPBS and F0 in formulating AAV.

Table A. Formulation Comparison

Component (unit)	DPBS	F0	F2	F3
NaCl (mM)	136.89	171.80	171.80	171.80
KCl (mM)	2.68	2.68	2.68	2.68
Na ₂ HPO ₄ (mM)	8.10	8.10	8.10	8.10
KH ₂ PO ₄ (mM)	1.47	1.47	1.47	1.47
CaCl ₂ (mM)	0.90	0.90	0.0	0.0
MgCl ₂ (mM)	0.50	0.50	0.49	1.40
Sucrose (% w/v)	0.0	1.00%	1.00%	1.00%
Poloxamer 188 (% w/v)	0.0	0.05%	0.05%	0.05%

[0050] As used herein, the concentrations of the various ingredients in the formulations may be expressed with zero, one, or two decimal places. Thus, for example, in F2, the concentrations of NaCl, KCl, Na₂HPO₄, KH₂PO₄, MgCl₂, and sucrose may be expressed respectively as 171.80 (or 171.8 or 172) mM, 2.68 (or 2.7) mM, 8.10 (or 8.1 or 8) mM, 1.47 (or 1.5) mM, 0.49 (or 0.5) mM, and 1.00% (or 1.0% or 1%) (w/v). In F3, the concentrations of NaCl, KCl, Na₂HPO₄,

KH₂PO₄, MgCl₂, and sucrose may be expressed respectively as 171.80 (or 171.8 or 172) mM, 2.68 (or 2.7) mM, 8.10 (or 8.1 or 8) mM, 1.47 (or 1.5) mM, 1.40 (or 1.4) mM, and 1.00% (or 1.0% or 1%) (w/v).

[0051] The pharmaceutical compositions may further comprise one or more preservatives such as ascorbic acid (vitamin C), sulfites, sorbates, benzoates, phenol, m-cresol, benzyl alcohol, benzalkonium chloride, phenoxyethanol, and/or parabens (e.g., methyl paraben). In some embodiments, the pharmaceutical compositions do not contain any added preservatives.

[0052] The pharmaceutical compositions may comprise also other reagents that enhance the effectiveness of the pharmaceutical composition. The pharmaceutical composition may contain delivery vehicles such as liposomes, nanocapsules, microparticles, microspheres, lipid particles, and vesicles.

III. Exemplary Recombinant AAV

[0053] In exemplary embodiments, the present disclosure provides improved pharmaceutical compositions comprising rAAV vector for factor VIII gene therapy. The rAAV vector, termed PF-07055480/SB-525 (or “SB-525” herein), is of the 2/6 pseudo serotype, which comprises an AAV6 capsid and a recombinant genome having AAV2 inverted terminal repeats (ITR). SB-525’s genome carries an expression cassette encoding a B-domain-deleted (BDD) version of human factor VIII (FVIII). *See, e.g.,* WO 2017/074526 (sequence #37). SB-525 is administered as an intravenous (IV) dose and possesses liver-specific tropism to provide long-term hepatic production of factor VIII protein in patients with hemophilia A. *See* WO 2020/028830. The secreted FVIII protein has the same amino acid sequence as approved recombinant anti-hemophilic factors (Refacto® and Xyntha®)

[0054] SB-525’s genome comprises an expression cassette for a human factor VIII BDD version, which has the amino acid sequence shown below:

```

MQIELSTCFFLCLLRFCFSATTRYYLGAVELSWDYMQSDLGELPVDARFPFPRVPKSFPF
NTSVVYKKTLEFVEFTDHLFNIAKPRPPWMGLLGPTIQAEVYDTVVITLKNMASHPVSLH
AVGVSYWKASEGAEYDDQTSQREKEDDKVFPGGSHTYVWQVLKENGPMASDPLCLTYSY
LSHVDLVKDLNSGLIGALLVCREGLAKEKTQTLHKFILLFAVFDEGKSWHSETKNSLM
QDRDAASARAWPKMHTVNGYVNRSLPGLIGCHRKSVYWHVIGMGTTPFVHSIFLEGHTF
LVRNHRQASLEISPITFLTAQTLMDLGQFLLFCHISSHQHDGMEAYVKVDSCPEEPQL
RMKNNEEAEDYDDDLTDSEMDVVRFDDDNSPSFIQIRSVAKKHPKTWVHYIAAEEDWD
YAPLV LAPDDRSYKSQYLNNGPQRIGRKYKKVRFMAYTDETFKTREAIQHESGILGPLL
YGEVGD TLLIIFKNQASRPYNIYPHGITDVRPLYSRRLPKG VKHLKDFPILPGEIFKYK

```

WTVTVEDGPTKSDPRCLTRYYSFVNMERDLASGLIGPLLCYKESVDQQRGNQIMSDKR
 NVILFSVFDENRSWYL TENIQRFLPNPAGVQLEDPEFQASNIMHSINGYVFDLSQLSVC
 LHEVAYWYILSIGAQTDFLSVFFSGYTFKHKMYEDTLTLFPFSGETVFMSENPGLWI
 LGCHNSDFRNRGMTALLKVSSCDKNTGDYYEDSYEDISAYLLSKNNAIEPRSFQONPPV
 LKRHQREITRTTLQSDQEEIDYDDTISVEMKKEDFDIYDEDENQSPRSFQKKTRHYFIA
 AVERLWDYGMSSSPHVLNRNAQSGSVPQFKKVVFQEFDTGDSFTQPLYRGELNEHLGLLG
 PYIRAEVEDNIMVTFRNQASRPYSFYSSLI SYEEDQRQGAEPKKNFVKPNETKTYFWKV
 QHHMAPTKDEFDCKAWAYFSDVDLEKDVHSGLIGPLLVCHTNTLNPAHGRQVTVQEFAL
 FFTIFDETKSWYFTENMERNCRAPCNIQMEDPTFKENYRFHAINGYIMDTLPGLVMAQD
 QRIRWYLLSMGSENENIHSIHFSGHVFTVRKKEEYKMALYNLYPGVFETVEMLPKAGI
 RVECLIGEHLHAGMSTLFLVYSNKCQTPLGMASGHIRDFQITASGQYGQWAPKLARLHY
 SGSINAWSTKEPFSWIKVDLLAPMIIHGKTQGARQKFSSLYISQFIIMYSLDGKKWQT
 YRGNSTGTLMVFFGNVDSSGIKHNI FNPPIIARYIRLHPHTHYSIRSTLRMELMGCDLNS
 CSMPPLGMESKAISDAQITASSYFTNMFATWSPSKARLHLQGRSNAWRPQVNNPKEWLQV
 DFQKTMKVTGVTTQGVKSLTSMYVKEFLISSQDGHQWTLFFQNGKVVFQGNQDSFT
 PVVNSLDPPLLTRYLRIHPQSWVHQIALRMEVLGCEAQDLY (SEQ ID NO:1)

[0055] The signal peptide portion of SEQ ID NO:1 is shown in box above and is cleaved off when the protein is secreted.

[0056] The SB-525 genome comprises the following nucleotide sequence:

1	CTGCGCGCTC	GCTCGCTCAC	TGAGGCCGCC	CGGGCAAAGC	CCGGGCGTCG	GGCGACCTTT
61	GGTCGCCCCG	CCTCAGTGAG	CGAGCGAGCG	CGCAGAGAGG	GAGTGGCCAA	CTCCATCACT
121	AGGGGTTCTT	GCGGCCTAAG	CTTGAACCA	TTGCCACCTT	CAGGGGGAGG	CTGCTGGTGA
181	ATATTAACCA	AGATCACCCC	AGTTACCGGA	GGAGCAAACA	GGGACTAAGT	TCACACGCGT
241	GGTACCGTCT	GTCTGCACAT	TTCGTAGAGC	GAGTGTTCCG	ATACTCTAAT	CTCCCTAGGC
301	AAGGTTTATA	TTTGTGTAGG	TTACTTATTC	TCCTTTTGTT	GAATAAGTCA	ATAATCAGAA
361	TCAGCAGGTT	TGGAGTCAGC	TTGGCAGGGA	TCAGCAGCCT	GGGTTGGAAG	GAGGGGGTAT
421	AAAAGCCCCCT	TCACCAGGAG	AAGCCGTCAC	ACAGATCCAC	AAGCTCCTGA	AGAGGTAAGG
481	GT'TTAAGTTA	TCGTTAGTTC	GTGCACCATT	AATGTTTAAT	TACCTGGAGC	ACCTGCCTGA
541	AATCATTTTT	TTTTTCAGGT	GGCTAGTATG	CAGATCGAGC	TCTCCACCTG	CTTCTTTCTG
601	TGCCTGTTGA	GATTCTGCTT	CAGCGCCACC	AGGAGATACT	ACCTGGGGGC	TGTGGAGCTG
661	AGCTGGGACT	ACATGCAGTC	TGACCTGGGG	GAGCTGCCTG	TGGATGCCAG	GTTCCCCCCC
721	AGAGTGCCCC	AGAGCTTCCC	CTTCAACACC	TCTGTGGTGT	ACAAGAAGAC	CCTGTTTGTG
781	GAGTTCACCTG	ACCACCTGTT	CAACATTGCC	AAGCCCAGGC	CCCCCTGGAT	GGGCCCTGCTG
841	GGCCCCACCA	TCCAGGCTGA	GGTGTATGAC	ACTGTGGTGA	TCACCCTGAA	GAACATGGCC
901	AGCCACCCTG	TGAGCCTGCA	TGCTGTGGGG	GTGAGCTACT	GGAAGGCCTC	TGAGGGGGCT
961	GAGTATGATG	ACCAGACCAG	CCAGAGGGAG	AAGGAGGATG	ACAAGGTGTT	CCCTGGGGGC
1021	AGCCACACCT	ATGTGTGGCA	GGTGTGAAG	GAGAATGGCC	CCATGGCCTC	TGACCCCTG
1081	TGCCTGACCT	ACAGCTACCT	GAGCCATGTG	GACCTGGTGA	AGGACCTGAA	CTCTGGCCTG
1141	ATTGGGGCCC	TGCTGGTGTG	CAGGGAGGGC	AGCCTGGCCA	AGGAGAAGAC	CCAGACCCTG
1201	CACAAGTTCA	TCCTGCTGTT	TGCTGTGTTT	GATGAGGGCA	AGAGCTGGCA	CTCTGAAACC
1261	AAGAACAGCC	TGATGCAGGA	CAGGGATGCT	GCCTCTGCCA	GGGCCTGGCC	CAAGATGCAC
1321	ACTGTGAATG	GCTATGTGAA	CAGGAGCCTG	CCTGGCCTGA	TTGGCTGCCA	CAGGAAGTCT
1381	GTGTACTGGC	ATGTGATTGG	CATGGGCACC	ACCCCTGAGG	TGCACAGCAT	CTTCTTGGAG
1441	GGCCACACCT	TCCTGGTCAG	GAACCACAGG	CAGGCCAGCC	TGGAGATCAG	CCCCATCACC
1501	TTCTGACTG	CCCAGACCCT	GCTGATGGAC	CTGGGCCAGT	TCCTGCTGTT	CTGCCACATC
1561	AGCAGCCACC	AGCATGATGG	CATGGAGGCC	TATGTGAAGG	TGGACAGCTG	CCCTGAGGAG
1621	CCCCAGCTGA	GGATGAAGAA	CAATGAGGAG	GCTGAGGACT	ATGATGATGA	CCTGACTGAC
1681	TCTGAGATGG	ATGTGGTGAG	GTTTGATGAT	GACAACAGCC	CCAGCTTCAT	CCAGATCAGG

1741	TCTGTGGCCA	AGAAGCACCC	CAAGACCTGG	GTGCACTACA	TTGCTGCTGA	GGAGGAGGAC
1801	TGGGACTATG	CCCCCTGGT	GCTGGCCCT	GATGACAGGA	GCTACAAGAG	CCAGTACCTG
1861	AACAATGGCC	CCCAGAGGAT	TGGCAGGAAG	TACAAGAAGG	TCAGGTTCAT	GGCCTACACT
1921	GATGAAACCT	TCAAGACCAG	GGAGGCCATC	CAGCATGAGT	CTGGCATCCT	GGGCCCCCTG
1981	CTGTATGGGG	AGGTGGGGGA	CACCCTGCTG	ATCATCTTCA	AGAACCAGGC	CAGCAGGCCC
2041	TACAACATCT	ACCCCCATGG	CATCACTGAT	GTGAGGCCCC	TGTACAGCAG	GAGGCTGCCC
2101	AAGGGGGTGA	AGCACCTGAA	GGACTTCCCC	ATCCTGCCTG	GGGAGATCTT	CAAGTACAAG
2161	TGGACTGTGA	CTGTGGAGGA	TGGCCCCACC	AAGTCTGACC	CCAGGTGCCT	GACCAGATAC
2221	TACAGCAGCT	TTGTGAACAT	GGAGAGGGAC	CTGGCCTCTG	GCCTGATTGG	CCCCCTGCTG
2281	ATCTGCTACA	AGGAGTCTGT	GGACCAGAGG	GGCAACCAGA	TCATGTCTGA	CAAGAGGAAT
2341	GTGATCCTGT	TCTCTGTGTT	TGATGAGAAC	AGGAGCTGGT	ACCTGACTGA	GAACATCCAG
2401	AGGTTCCTGC	CCAACCCCTGC	TGGGGTGCAG	CTGGAGGACC	CTGAGTTCCA	GGCCAGCAAC
2461	ATCATGCACA	GCATCAATGG	CTATGTGTTT	GACAGCCTGC	AGCTGTCTGT	GTGCCTGCAT
2521	GAGGTGGCCT	ACTGGTACAT	CCTGAGCATT	GGGGCCAGA	CTGACTTCCT	GTCTGTGTTC
2581	TTCTCTGGCT	ACACCTTCAA	GCACAAGATG	GTGTATGAGG	ACACCCTGAC	CCTGTTCCCC
2641	TTCTCTGGGG	AGACTGTGTT	CATGAGCATG	GAGAACCCTG	GCCTGTGGAT	TCTGGGCTGC
2701	CACAACCTCTG	ACTTCAGGAA	CAGGGGCATG	ACTGCCCTGC	TGAAAGTCTC	CAGCTGTGAC
2761	AAGAACACTG	GGGACTACTA	TGAGGACAGC	TATGAGGACA	TCTCTGCCTA	CCTGCTGAGC
2821	AAGAACAATG	CCATTGAGCC	CAGGAGCTTC	AGCCAGAATC	CACCCGTCCT	TAAGCGCCAT
2881	CAGCGCGAGA	TCACCAGGAC	CACCCTGCAG	TCTGACCAGG	AGGAGATTGA	CTATGATGAC
2941	ACCATCTCTG	TGGAGATGAA	GAAGGAGGAC	TTTGACATCT	ACGACGAGGA	CGAGAACCAG
3001	AGCCCCAGGA	GCTTCCAGAA	GAAGACCAGG	CACTACTTCA	TTGCTGCTGT	GGAGAGGCTG
3061	TGGGACTATG	GCATGAGCAG	CAGCCCCCAT	GTGCTGAGGA	ACAGGGCCCA	GTCTGGCTCT
3121	GTGCCCCAGT	TCAAGAAGGT	GGTGTTCAG	GAGTTCACTG	ATGGCAGCTT	CACCCAGCCC
3181	CTGTACAGAG	GGGAGCTGAA	TGAGCACCTG	GGCCTGCTGG	GCCCCTACAT	CAGGGCTGAG
3241	GTGGAGGACA	ACATCATGGT	GACCTTCAGG	AACCAGGCCA	GCAGGCCCTA	CAGCTTCTAC
3301	AGCAGCCTGA	TCAGCTATGA	GGAGGACCAG	AGGCAGGGGG	CTGAGCCCAG	GAAGAACTTT
3361	GTGAAGCCCA	ATGAAACCAA	GACCTACTTC	TGGAAGGTGC	AGCACCACAT	GGCCCCCACC
3421	AAGGATGAGT	TTGACTGCAA	GGCCTGGGCC	TACTTCTCTG	ATGTGGACAT	GGAGAAGGAT
3481	GTGCACTCTG	GCCTGATTGG	CCCCCTGCTG	GTGTGCCACA	CCAACACCCT	GAAACCCTGCC
3541	CATGGCAGGC	AGGTGACTGT	GCAGGAGTTT	GCCCTGTTCT	TCACCATCTT	TGATGAAACC
3601	AAGAGCTGGT	ACTTCACTGA	GAACATGGAG	AGGAAGTGA	GGGCCCCCTG	CAACATCCAG
3661	ATGGAGGACC	CCACCTTCAA	GGAGAAGTAC	AGGTTCCATG	CCATCAATGG	CTACATCATG
3721	GACACCCTGC	CTGGCCTGGT	GATGGCCAG	GACCAGAGGA	TCAGGTGGTA	CCTGCTGAGC
3781	ATGGGCAGCA	ATGAGAACAT	CCACAGCATC	CACCTTCTCTG	GCCATGTGTT	CACTGTGAGG
3841	AAGAAGGAGG	AGTACAAGAT	GGCCCTGTAC	AACCTGTACC	CTGGGGTGTG	TGAGACTGTG
3901	GAGATGCTGC	CCAGCAAGGC	TGGCATCTGG	AGGGTGGAGT	GCCTGATTGG	GGAGCACCTG
3961	CATGCTGGCA	TGAGCACCCCT	GTTCTCTGGT	TACAGCAACA	AGTGCCAGAC	CCCCCTGGGC
4021	ATGGCCTCTG	GCCACATCAG	GGACTTCCAG	ATCACTGCCT	CTGGCCAGTA	TGGCCAGTGG
4081	GCCCCCAAGC	TGGCCAGGCT	GCACTACTCT	GGCAGCATCA	ATGCCTGGAG	CACCAAGGAG
4141	CCCTTCAGCT	GGATCAAGGT	GGACCTGCTG	GCCCCCATGA	TCATCCATGG	CATCAAGACC
4201	CAGGGGGCCA	GGCAGAAGTT	CAGCAGCCTG	TACATCAGCC	AGTTCATCAT	CATGTACAGC
4261	CTGGATGGCA	AGAAGTGGCA	GACCTACAGG	GGCAACAGCA	CTGGCACCCCT	GATGGTGTTC
4321	TTTGGCAATG	TGGACAGCTC	TGGCATCAAG	CACAACATCT	TCAACCCCCC	CATCATTGCC
4381	AGATACATCA	GGCTGCACCC	CACCCACTAC	AGCATCAGGA	GCACCCTGAG	GATGGAGCTG
4441	ATGGGCTGTG	ACCTGAACAG	CTGCAGCATG	CCCCTGGGCA	TGGAGAGCAA	GGCCATCTCT
4501	GATGCCCAGA	TCACTGCCAG	CAGCTACTTC	ACCAACATGT	TTGCCACCTG	GAGCCCCAGC
4561	AAGGCCAGGC	TGCATCTGCA	GGGCAGGAGC	AATGCCTGGA	GGCCCCAGGT	CAACAACCCC
4621	AAGGAGTGGC	TGCAGGTGGA	CTTCCAGAAG	ACCATGAAGG	TGACTGGGGT	GACCACCCAG
4681	GGGGTGAAGA	GCCTGCTGAC	CAGCATGTAT	GTGAAGGAGT	TCCTGATCAG	CAGCAGCCAG
4741	GATGGCCACC	AGTGGACCCT	GTTCTTCCAG	AATGGCAAGG	TGAAGGTGTT	CCAGGGCAAC
4801	CAGGACAGCT	TCACCCCTGT	GGTGAACAGC	CTGGACCCCC	CCCTGCTGAC	CAGATACCTG
4861	AGGATTACCC	CCCAGAGCTG	GGTGCACCAG	ATTGCCCTGA	GGATGGAGGT	GCTGGGCTGT
4921	GAGGCCCAGG	ACCTGTACTG	AGGATCCAAT	AAAATATCTT	TATTTTCATT	ACATCTGTGT
4981	GTTGGTTTTT	TGTGTGTTTT	CCTGTAACGA	TCGGGCTCGA	GCGCAGGAAC	CCCTAGTGAT
5041	GGAGTTGGCC	ACTCCCTCTC	TGCGCGCTCG	CTCGCTCACT	GAGGCCGCC	GGGCTTTGCC

5101 CGGGCGGCCT CAGTGAGCGA GCGAGCGCGC AG (SEQ ID NO:2)

[0057] In the above sequence, the left (5') ITR (AAV2 ITR) spans nucleotides 1-130, and the right (3') ITR (AAV2 ITR) spans nucleotides 5025-5132. Both ITRs are boxed.

IV. Use of Recombinant AAV Formulations

[0058] The pharmaceutical compositions of the invention may be supplied in an article of manufacture (e.g., a kit) comprising vials (e.g., pre-treated glass vials or COP vials) and instructions for use. In some embodiments, each vial contains about 1E+11 to 1E+15 vg per mL of AAV in 0.5-50 mL (e.g., 1-10 mL). In particular embodiments, each vial contains 5E+12 to 1E+14/mL (e.g., 1E+13 vg/mL). In certain embodiments, each vial contains 6E+13 vg in 6 mL.

[0059] The compositions may be administered to the patients once or more than once. For example, the compositions may be administered to the patients with an interval of no less than 1 month, 3 months, 6 months, 9 months, or one year. In some embodiments, the composition may be administered to the patients with an interval of no less than two, five, seven, ten, or fifteen years. The pharmaceutical composition may be provided to a patient in need thereof through a route appropriate for the disease to be treated. For example, the compositions may be administered through intravenous injection, intraarterially injection, intracranial injection, intraperitoneal injection, portal vein injection, or intramuscular injection. For example, the SB-525 pharmaceutical compositions may be provided intravenously to hemophilia A patients at 1E+11 to 1E+15 vg/kg, such as 1E+11 to 1E+14 (e.g., 1E+12 to 1E+14) vg/kg. In some embodiments, the SB-525 pharmaceutical compositions may be provided intravenously to hemophilia A patients at 5E+11, 6E+11, 7E+11, 8E+11, 9E+11, 1E+12, 2E+12, 3E+12, 4E+12, 5E+12, 6E+12, 7E+12, 8E+12, 9E+12, 1E+13, 2E+13, 3E+13, 4E+13, 5E+13, 6E+13, 7E+13, 8E+13, 9E+13, or 1E+14 vg/kg. In some embodiments, the SB-525 compositions are provided intravenously to hemophilia A patients at a dose of about 6E+13 vg/kg.

[0060] In some embodiments, the patients have with severe or moderate hemophilia A. In further some embodiments, the patients have no inhibitors (alloantibodies to Factor VIII). In certain embodiments, the patients do not have neutralizing antibodies to AAV6. The patients may be adult or adolescent patients (≥ 12 years of age) or pediatric patients (< 12 years of age).

[0061] Unless otherwise defined herein, scientific and technical terms used in connection with the present disclosure shall have the meanings that are commonly understood by those of ordinary skill in the art. Exemplary methods and materials are described below, although methods and materials similar or equivalent to those described herein can also be used in the practice or testing of the present disclosure. In case of conflict, the present specification, including definitions, will control. Generally, nomenclature used in connection with, and techniques of, cardiology, medicine, medicinal and pharmaceutical chemistry, and cell biology described herein are those well-known and commonly used in the art. Enzymatic reactions and purification techniques are performed according to manufacturer's specifications, as commonly accomplished in the art or as described herein. Further, unless otherwise required by context, singular terms shall include pluralities and plural terms shall include the singular. Throughout this specification and embodiments, the words "have" and "comprise," or variations such as "has," "having," "comprises," or "comprising," will be understood to imply the inclusion of a stated integer or group of integers but not the exclusion of any other integer or group of integers. All publications and other references mentioned herein are incorporated by reference in their entirety. Although a number of documents are cited herein, this citation does not constitute an admission that any of these documents forms part of the common general knowledge in the art. As used herein, the term "approximately" or "about" as applied to one or more values of interest refers to a value that is similar to a stated reference value. In certain embodiments, the term refers to a range of values that fall within 10%, 9%, 8%, 7%, 6%, 5%, 4%, 3%, 2%, 1%, or less in either direction (greater than or less than) of the stated reference value unless otherwise stated or otherwise evident from the context.

[0062] In order that this invention may be better understood, the following examples are set forth. These examples are for purposes of illustration only and are not to be construed as limiting the scope of the invention in any manner.

EXAMPLES

[0063] The following Examples describe a study in which the inventors discovered a current SB-525 formulation formed precipitates upon freeze/thaw cycling, a common test for assessing the shelf-life and long-term stability of biological products. The inventors discovered that precipitates also formed in the buffer formulation alone (i.e., buffer formulation containing no

active viral vector ingredient). The inventors thus set out to find a new formulation that would circumvent this problem.

[0064] In the study, the SB-525 viral preparations were prepared and formulated at approximately $1.0\text{E}+13$ vg/mL in phosphate-buffered saline (PBS) with the addition of CaCl_2 , MgCl_2 , 35 mM NaCl, 1% sucrose, and 0.05% Kolliphor® P 188 (poloxamer 188), filled at 5 mL into 6 mL Aseptic Technologies crystal closed vials, and stored at $\leq -65^\circ\text{C}$. For the stability tests, SB-525 compositions were placed on stability at -0°C , 5°C , and 25°C for up to 24 months (**Table 5**). Five cycles of uncontrolled freeze/thaw (F/T) (**Table 1**) and a 24-hour agitation (AG) (**Table 2**) were also performed. Drug product (DP) quality attributes were analyzed according to **Table 5** to **Table 8**.

Example 1: Freeze/Thaw Cycling and Agitation Tests

Materials and Methods

[0065] The SB-525 was provided at $1.0\text{E}+13$ vg/mL nominal in a formulation buffer containing the following ingredients: 0.90 mM CaCl_2 , 0.49 mM MgCl_2 , 2.68 mM KCl, 1.47 mM KH_2PO_4 , 172 mM NaCl, 8.10 mM Na_2HPO_4 , 1% (w/v) sucrose, 0.05% (w/v) poloxamer 188, pH 7.36. 23 vials were stored at -70°C , 11 at 5°C , and 8 at 25°C .

[0066] For freeze/thaw (F/T) cycling tests, samples were uncontrolled F/T cycled for three and five cycles as shown in **Table 1**.

Table 1. DP Uncontrolled F/T Cycling

F/T Cycling	Day 1	Day 2		Day 3		Day 4		Day 5	
	4 pm	8 am	4 pm	8 am	4 pm	8 am	4 pm	8 am	4 pm
3	-80°C	25°C	-80°C	25°C	-80°C	25°C			
5	-80°C	25°C	-80°C	25°C	-80°C	25°C	-80°C	25°C	-80°C

[0067] For agitation tests, samples were agitated at room temperature as shown in **Table 2**.

Table 2. DP Agitation Setup

Agitation Stability	T_0	$T_{24\text{Hr}}$
	Day 1 (8 am)	Day 2 (8 am)
300 rpm	Use Stability T_0	Pull 24Hr material

[0068] The vials stored at -70°C, 5°C, at 25°C, respectively, were analyzed for (i) appearance (liquid), pH, dynamic light scattering (DLS), and high accuracy liquid particle counter (HIAC); (ii) vg identity (qPCR), vg titer (qPCR), capsid titer (ELISA), capsid identity (ELISA), and infectious titer (TCID₅₀); (iii) *in vitro* FVIII Activity (bioassay); or (iv) UV 260/280, reduced CE-SDS, SEC titer_{260/280}.

[0069] Samples that went through F/T cycling were analyzed for (i) appearance (liquid), pH, osmolarity, DLS, HIAC; (ii) vg titer (qPCR), capsid titer and identity (ELISA), reduced SDS-PAGE, and infectious titer (TCID₅₀); or (iii) UV 260/280, SEC titer_{260/280}, RP-HPLC, *in vitro* FVIII Activity (bioassay), and reduced CE-SDS.

[0070] Samples that went through 24-hour agitation (AG) were analyzed for (i) appearance (liquid), pH, osmolarity, DLS, and HIAC; (ii) vg titer (qPCR), capsid titer and identity (ELISA), reduced SDS-PAGE, infectious titer (TCID₅₀); or (iii) UV 260/280, SEC titer_{260/280}, RP-HPLC, *in vitro* FVIII activity (bioassay), and reduced CE-SDS.

Results

[0071] Results from the SB-525 DP stability, 5 cycles of uncontrolled F/T, and 24-hour agitation arms are described below.

[0072] As shown in **Table 3** below, there was no significant change in color, clarity, pH, hydrodynamic radius as measured by DLS, vg titer, capsid titer, capsid purity as measured by SDS-PAGE, infectious titer, UV 260/280 (a surrogate measure for empty:full capsid ratio), and biological activity (% relative potency) when the samples were stored at the tested conditions for up to one month. There was a change in visible particles at 5°C and 25°C. All the samples were frozen and appeared white and opaque at T zero (T₀).

Table 3. Stability Results

Test		Storage	T ₀ [*]	T _{1M} ²
Appearance, Liquid	Color	-70°C	Colorless, clear, no precipitates, no contaminants	B6
		5°C		B6
		25°C		B6
	Clarity	-70°C		<= Ref. I
		5°C		<= Ref. I
		25°C		<= Ref. I
	Visible Particles	-70°C		EFVP
		5°C		1-5 particles
		25°C		TMTC
pH		-70°C	7.3	7.3

Test		Storage	T ₀ *	T _{1M} ²
		5°C		7.3
		25°C		7.1
		-70°C		1.70e ⁻⁷
DLS	Z-average – Diffusion Coefficient (cm ² /s)	5°C	-	1.70e ⁻⁷
		25°C		1.70e ⁻⁷
		-70°C		1.70e ⁻⁷
	Z-average – Radius (nm)	-70°C	22.3 (diameter)	14.7
		5°C		14.8
		25°C		14.4
	Peak 1 – Diffusion Coefficient (cm ² /s)	-70°C	-	1.74e ⁻⁷
		5°C		1.69e ⁻⁷
		25°C		1.70e ⁻⁷
	Peak 1 – radius (nm)	-70°C		14.2
		5°C		14.6
		25°C		14.5
	Peak 2 – Diffusion Coefficient (cm ² /s)	-70°C		9.42e ⁻¹⁰
		5°C		1.26e ⁻⁹
		25°C		-
	Peak 2 – radius (nm)	-70°C		2658.7
		5°C		2261.8
		25°C		-
Vg titer by qPCR (vg/mL) ¹		-70°C	8.74E+12	9.77E+12
		5°C		1.01E+13
		25°C		9.98E+12
Capsid ELISA (capsids/mL)		-70°C	2.76E+13	2.52E+13
		5°C		2.53E+13
		25°C		2.66E+13
SDS-PAGE		-70°C	VP1, VP2, and VP3 comparable to reference. No other bands ≥ sensitivity control observed (% Purity VP1, VP2, VP3 combined) = 100%	VP1, VP2, and VP3 comparable to reference. No other bands ≥ sensitivity control observed (% Purity VP1, VP2, VP3 combined) = 100%
		5°C		
		25°C		
Infectious Titer TCID ₅₀ (TCID ₅₀ /mL)		-70°C	1.12E+9	3.56E+9
		5°C		2.00E+9
		25°C		1.12E+9
SEC UV 260/280 Ratio		-70°C	1.37	1.36
		5°C		1.36
		25°C		1.36
<i>In vitro</i> FVIII Activity/ Functional Bioassay (% Relative Potency)		-70°C	85 ³	80
		5°C		80
		25°C		69

*T₀: all analytics, except functional bioassay, performed per release testing

¹ T₀ re-test of Vector Genome Identity & Titer by qPCR was performed

² Study terminated after T1M due to particle detection in DP; investigation was initiated

³ T₀ re-test for Bioassay (% Relative Potency) against a different reference material

EFVP: Essentially Free of Visible Particles

TMTC: Too many to Count

N.T.: Not tested

[0073] As shown in **Table 4**, there was no significant change in capsid purity after one month at all storage conditions.

Table 4. Capsid Purity Stability Results

Time Point	Storage	Capsid (%)	VP1 (%)	VP2 (%)	VP3 (%)	Other (%)
T ₀		94.3	2.2	3.4	94.4	5.7
T _{1M}	-70°C	95.1	1.1	3.1	95.8	4.9
	5°C	95.0	2.1	3.2	94.7	5.0
	25°C	94.8	1.9	3.0	95.1	5.2

[0074] As shown in **Table 5**, there was no significant aggregation or formation of high molecular weight species (HMWS) after one month at all storage conditions (NMT: not more than; LOQ: limit of quantification).

Table 5. SEC Titer Stability Results After One Month
Under Intended, Accelerated, and Stressed Storage Conditions

Time Point	Storage	Viral Particle (vp) Titer (vp/mL)	Monomer (%)	HMWS (%)	LOQ
T ₀		5.12E+13	100	NMT 1.0	1.01E+12
T _{1M}	-70°C	2.65E+13	100	NMT 1.0	1.01E+12
	5°C	2.57E+13	100	NMT 1.0	1.01E+12
	25°C	3.31E+13	100	NMT 1.0	1.01E+12

[0075] As shown in **Table 6**, there was no significant increase in subvisible particles (per HIAC) after one month at all storage conditions. All time points showed particles below the USP <787> standard in regard to the $\geq 10 \mu\text{m}$ and the $\geq 25 \mu\text{m}$ ranges, where acceptance criteria are no more than 6,000 particles $> 10 \mu\text{m}$ per container and no more than 600 particles $> 25 \mu\text{m}$ per container.

Table 6. HIAC Stability Results

Storage	T ₀ (particles/container)		T _{1M} (particles/mL)				
	$\geq 10 \mu\text{m}$	$\geq 25 \mu\text{m}$	$\geq 2 \mu\text{m}$	$\geq 5 \mu\text{m}$	$\geq 8 \mu\text{m}$	$\geq 10 \mu\text{m}$	$\geq 25 \mu\text{m}$
-70°C	10	3	143	11	1	0	0
5°C			626	96	24	15	0
25°C			478	94	15	4	1

[0076] As shown in **Table 7**, there was no significant change in color, clarity, pH, hydrodynamic radius as measured by DLS, vg titer, capsid titer, infectious titer, UV 260/280 (surrogate for empty:full capsid ratio), capsid purity as measured by RP-HPLC, and biological activity (% relative potency) upon F/T/AG stress. There was a change in visible particles for F/T.

Table 7. Cycle Uncontrolled F/T & 24-Hour AG Results

Test		Treatment	Results
Appearance	Color	F/T	B9 ¹
		AG	B9
	Clarity	F/T	<= REF 1 ²
		AG	<= REF 1
	Visible Particles	F/T	TMTC
		AG	Essentially free of visible particles (EFVP)
pH		F/T	7.2
		AG	7.3
Osmolality (mOsm/kg)		F/T	387
		AG	388
DLS	Z-average – Diffusion Coefficient (cm ² /s)	F/T	1.70e ⁻⁷
		AG	1.67e ⁻⁷
	Z-average –Radius (nm)	F/T	14.5
		AG	14.7
	Peak 1 – Diffusion Coefficient (cm ² /s)	F/T	1.69e ⁻⁷
		AG	1.56e ⁻⁷
	Peak 1 – radius (nm)	F/T	14.5
		AG	15.7
Vg titer by qPCR (vg/mL)		F/T	8.44E+12
		AG	8.63E+12
Capsid ELISA (capsids/mL)		F/T	2.08E+13
		AG	2.00E+13
SDS-PAGE		F/T	VP1, VP2, and VP3 comparable to reference. No other bands ≥ sensitivity control observed (% Purity VP1, VP2, VP3 combined) = 100%
		AG	
Infectious Titer TCID50 (TCID ₅₀ /mL)		F/T	9.29E+8
		AG	2.42E+9
UV 260/280 Ratio		F/T	1.32
		AG	1.32
SEC UV 260/260 Ratio		F/T	1.35
		AG	1.35
RP-HPLC	Capsid (%)	F/T	98.9

Test	Treatment	Results
	AG	99.0
	F/T	1.1
	AG	1.0
<i>in vitro</i> FVIII Activity/ Functional Bioassay (% Relative Potency)	F/T	85
	AG	80

¹ Color is not more intense or more opalescent than reference standard listed.

² Clarity is not more opalescent than reference standard listed.

[0077] As shown in **Table 8**, there was no significant change in reduced CE-SDS upon F/T/AG stress.

Table 8. 5 Cycle Uncontrolled F/T & 24-Hour AG CE-SDS-Reduced Results

Treatment	Capsid (%)	VP1 (%)	VP2 (%)	VP3 (%)	Other (%)
F/T	95.1	2.0	3.0	95.0	4.9
AG	95.2	2.0	3.2	94.8	4.8

[0078] As shown in **Table 9**, there was no significant aggregation/HMWS formation upon F/T/AG stress.

Table 9. 5 Cycle Uncontrolled F/T & 24-Hour AG SEC Titer Results

Treatment	Viral Particle Titer (vp/mL)	Monomer (%)	HMWS (%)	LOQ
F/T	4.95E+13	100	0	1.01E+12
AG	5.05E+13	100	0	1.01E+12

[0079] As shown in **Table 10**, there was a significant increase in subvisible particles (per HIAC) upon F/T. Particles upon F/T were above USP <787> standard with regard to the $\geq 10 \mu\text{m}$ and the $\geq 25 \mu\text{m}$ ranges; whereas particles upon AG were below USP <787> standard with regard to the $\geq 10 \mu\text{m}$ and the $\geq 25 \mu\text{m}$ ranges.

Table 10. 5 Cycle Uncontrolled F/T & 24-Hour AG HIAC Results

Treatment	$\geq 2 \mu\text{m}$ (particles/mL)	$\geq 5 \mu\text{m}$ (particles/mL)	$\geq 8 \mu\text{m}$ (particles/mL)	$\geq 10 \mu\text{m}$ (particles/mL)	$\geq 25 \mu\text{m}$ (particles/mL)
F/T	24551	14051	8725	6576	710
AG	439	41	12	4	0

[0080] The above results demonstrate that the SB-525 samples showed no significant change in color, clarity, pH, hydrodynamic radius (as measured by dynamic light scattering), vg titer, capsid titer, infectious titer, empty:full capsid ratio (as measured by UV 260/280), biological activity (% relative potency), capsid purity (as measured by SDS-PAGE, reduced CE-SDS and RP-HPLC), or viral particle titer when stored up to one month at -70°C, 5°C, or 25°C.

Additionally, there was no significant change in color, clarity, pH, hydrodynamic radius (as measured by dynamic light scattering), vg titer, capsid titer, infectious titer, empty:full capsid ratio (as measured by UV 260/280), biological activity (% relative potency), capsid purity (as measured by SDS-PAGE, reduced CE-SDS, RP-HPLC), or viral particle titer upon F/T or AG stress, either. However, there was a significant increase in subvisible particles (per HIAC) upon F/T and an increase in visible particles after one month of storage at 5°C and 25°C. The study was terminated at one month because of these results.

Example 2: Characterization of Particulate Matter in DP Vials after F/T

[0081] The samples from the above stability tests of SB-525 formulations, in which particulate matter was observed, were analyzed further for isolation and identification.

[0082] Drug product samples that had undergone 5 freeze-thaw cycles (“DP after 5 F/T”) and samples that had been stored at 25°C for one month (DP after 1M at 25°C) were analyzed. Particulate matter was isolated onto 0.8 µm gold filters. It was then imaged on the filter under partial ring lighting using a Keyence VHX6000 digital microscope at 150x magnification. The filters were then transferred to the Fourier transform infrared (FTIR) microscope where spectra were obtained for the material. FTIR spectra were compared against the KnowItAll database for known spectra. A portion of the DP after 1M at 25°C sample was also scraped off the filter onto a glass slide and imaged under plane polarized light and crossed polars using a Nikon Eclipse ME600 polarized light microscope. A portion of the gold filter was cut and analyzed under the JEOL6000 SEM equipped with an EDS module. The sample was also analyzed by Raman microscopy analysis.

[0083] For all samples, digital imaging showed a white solid that appeared to be semi-crystalline. Under crossed polars, the “DP after 5 F/T” sample did not exhibit any birefringence, suggesting it was either an amorphous substance or isotropic. SEM/EDS analysis showed large portions of oxygen, phosphorous, and calcium. FTIR analysis identified the material as a salt,

but further identification was required. Raman microscopy was able to detect peaks characteristic of calcium phosphate, suggesting the white particulate matter in all samples was calcium phosphate. **FIGs. 1A-F** summarize the particle characterization studies.

Example 3: Reformulation of SB-525

[0084] This example describes experiments testing new SB-525 formulations that can circumvent the precipitation problem observed above. These experiments assessed the short-term stability of SB-525 DP when reformulated to (1) remove calcium, (2) remove calcium and magnesium, or (3) increase sucrose to 8.5%. These changes were expected to prevent generation of particles by (1) removing the source of particles, (2) removing both divalent cations due to solubility concerns of MgCl_2 , or (3) increasing stability of the formulation against freeze-thaw stress, respectively. Two pilot batches were buffer exchanged and filled at 2.5 mL into 6 mL AT vials to mimic worst-case conditions for surface area to volume (SA/V) fill. After fill/finish, vials were subjected to 5 uncontrolled freeze-thaw cycles ($\leq -65^\circ\text{C}$ to ambient), and then placed on stability at $\leq -65^\circ\text{C}$ (intended storage), $2-8^\circ\text{C}$ (stressed; liquid storage), 25°C (accelerated; liquid storage), and 40°C (aggressive/forced degradation condition; liquid storage). These conditions were selected with the expectation of highlighting differences between formulations.

[0085] SB-525 virus was purified from the clarified bulk harvest and formulated at approximately $1.0\text{E}+13$ vg/mL in phosphate buffered saline (PBS) containing CaCl_2 , MgCl_2 , 35 mM NaCl, 1% Sucrose, and 0.05% Kolliphor® P 188 (poloxamer 188), filled into vials, and stored at $\leq -65^\circ\text{C}$.

Materials and Methods

[0086] The SB-525 formulated bulk drug substance contained: SB-525 ($1.0\text{E}+13$ vg/ mL nominal) in 0.90 mM CaCl_2 , 0.49 mM MgCl_2 , 2.68 mM KCl, 1.47 mM KH_2PO_4 , 172 mM NaCl, 8.10 mM Na_2HPO_4 , 1% (w/v) sucrose, 0.05% (w/v) poloxamer 188, pH 7.0-7.6, sterile filtered and stored in 125 mL HDPE bottles at $\leq -65^\circ\text{C}$. Na_2HPO_4 (sodium phosphate, dibasic, anhydrous), KH_2PO_4 (potassium phosphate monobasic), CaCl_2 (calcium chloride dihydrate), and MgCl_2 (magnesium chloride hexahydrate) were sourced from JT Baker or Fischer. Poloxamer 188 was sourced from BASF. Storage vials were 6 mL Aseptic Technologies (AT) Closed Crystal Vials (Aseptic Technologies; cat# VIA-060000), which are vials with primary container closure, cyclo-olefin copolymer (COC) with in-place thermoplastic elastomer stopper and yellow

cap. Prior to DP fill, the formulated bulk drug substance (FBDS) was filled into 125 mL high density polyethylene (HDPE) bottles and frozen to $\leq -65^{\circ}\text{C}$. Bottles were then thawed at ambient temperature, buffer exchanged, filtered, and DP filled into 6 mL AT vials using the Aseptic Technologies M1 unit. One pilot batch material was buffer exchanged via TFF using 50 kDa filter units at a pressure of approximately 40 psi for a total of 10 exchange volumes. The second pilot batch was buffer exchanged using Amicon stir cells over a 50 kDa NMW PES filter at approximately 40 psi for a total of 10 exchange volumes. The two exchange methods have been previously used for AAV without significant loss or adsorption of material.

[0087] The AT vials were filled 2.5 mL at a target of approximately $1.0\text{E}+13$ vg/mL, freeze thaw cycled 5 times at an uncontrolled rate from $\leq -65^{\circ}\text{C}$ to ambient conditions, and placed on stability in non-GMP storage units. Vials were subjected to the stability pull schedule outlined in **Table 11**, and tested by the methods listed in **Table 12**. Formulations tested are outlined in **Table 13**.

Table 11. Capsid Purity Stability Schedule

Storage Temp ($^{\circ}\text{C}$)	T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}	T _{5W}
-70 (± 10)	XABCD	-	-	-	-	XABC	XA
5 (± 3)		X	X	XA	X		XA
25 (± 3)		X	X	XAB	X		-
40 (± 3)		XABCF	XA	X	XF	-	-

Table 12. Analytical Testing

Test Code	Test Method
X	Appearance (Liquid)
A	pH
A	SEC – Aggregates, UV 260/280, Viral Particle Titer
B	Functional FVIII Bioassay
C	Light Obscuration with HIAC
D	Osmolality (freezing point depression)
D	P188
D	Sucrose
D	Reduced capillary gel electrophoresis (rCGE) (ratio and purity)
F	Mass Spec (deamidation)
B	Vector Genome Titer (qPCR)

Table 13. Formulations

Exchange Process	Material	Formulation	Difference from Control
	M05	Control	N/A

Exchange Process	Material	Formulation	Difference from Control
Millipore Amicon Stir Cell	M06	Ca Negative	Removed CaCl ₂
	M07	Ca /Mg Negative	Removed CaCl ₂ and MgCl ₂
	M08	8.5% Sucrose	Increase sucrose to 8.5% (w/v)

*Control formulation (F0): 0.90 mM CaCl₂, 0.49 mM MgCl₂, 2.68 mM KCl, 1.47 mM KH₂PO₄, 172 mM NaCl, 8.10 mM Na₂HPO₄, 1% (w/v) sucrose, 0.05% (w/v) poloxamer 188 (P188), pH 7.0-7.6.

Results

Formulation Characterization

[0088] To confirm that the correct formulations were achieved, SB-525 DP formulations were tested to determine osmolality, P188 concentration, sucrose concentration, and capsid ratio and purity by reduced CGE. The results to these tests for M05-M08 are shown in **Table 14** below. Capsid purity and ratio showed no significant differences as measured by rCGE. Osmolality is within expected ranges, dependent on the level of sucrose. M05-M08 samples demonstrated the correct concentrations for P188. In addition, concentrations of calcium and magnesium were confirmed by CEDEX to be within expected ranges for all formulations.

Table 14. Formulation Characterization

Material	Osmolality (mOsm/kg)	Sucrose (%)	P188 (µg/mL)	Purity (%Capsid)	AAV Ratio		
					%VP1	%VP2	%VP3
M05 – Control	399	1.10	589	95.2	2.0	3.3	94.7
M06 –Ca Negative	391	1.00	565	95.3	2.1	3.0	94.9
M07 – Ca/Mg Negative	417	1.00	575	95.8	1.7	2.9	95.4
M08 – 8.5% Sucrose	655	8.50	568	95.2	1.8	2.6	95.6

Vector Genome Titer (vg/mL)

[0089] Vector genome titer results for M05-M08 are presented below in **Table 15**. At $\leq -65^{\circ}\text{C}$ and $2-8^{\circ}\text{C}$ conditions, there was no significant trend observed after storage for 4 weeks. After 40°C for 3 days, the Ca/Mg negative formulation showed a downward trend in vector genome titer. This result is considered within assay variability, but is consistent with vp titer and UV 260/280 results for this formulation.

Table 15. Genome Titer (vg/mL)

Formulation	Timepoint					
	T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}
-70°C						
M05 - Control	7.44E+12	Not tested				7.77E+12

Formulation	Timepoint					
	T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}
M06 - Ca Negative	7.98E+12					7.47E+12
M07 - Ca / Mg Negative	6.33E+12					6.24E+12
M08 - 8.5% Sucrose	6.23E+12					6.56E+12
5°C						
M05 - Control	7.44E+12	Not tested				7.39E+12
M06 - Ca Negative	7.98E+12					6.95E+12
M07 - Ca / Mg Negative	6.33E+12					6.56E+12
M08 - 8.5% Sucrose	6.23E+12					6.24E+12
25°C						
M05 - Control	7.44E+12	Not tested		7.23E+12	Not tested	7.84E+12
M06 - Ca Negative	7.98E+12			7.41E+12		7.09E+12
M07 - Ca / Mg Negative	6.33E+12			6.04E+12		5.85E+12
M08 - 8.5% Sucrose	6.23E+12			6.11E+12		7.01E+12
40°C						
M05 - Control	7.44E+12	7.16E+12	Not tested			
M06 - Ca Negative	7.98E+12	6.94E+12				
M07 - Ca / Mg Negative	6.33E+12	4.63E+12				
M08 - 8.5% Sucrose	6.23E+12	6.68E+12				

Visible Appearance

[0090] All formulations displayed B9 color and less than or equal to opalescence reference 1, color and clarity standards, respectively, used during visual inspection against a black and white background, for all timepoints and conditions (*see, e.g., Ph. Eur. 7.0, 20201, 20202 (01/2008); Millipore Sigma Colour Reference Solutions B*). Visible appearance results show that only the “control” formulation (M05) displayed white flake-like particles that were too many to count (TMTC) when held at 25°C for longer than 1 week, or 40°C for 3 days or longer. M06 and M07 formulations did occasionally show TMTC visible particles; however, the particle descriptions for these samples were all fiber-like particles. While these particles were not identified for this study, fiber-like particles are more characteristic of exogenous particles (*e.g., filter particles*) versus the flake-like particles that are characteristic of calcium phosphate. Also, appearance of fiber-like/exogenous particles is negligible for lab-based development studies where DP is not manufactured under tightly controlled environmental conditions.

Cumulative Sub-Visible Particles

[0091] Sub-visible cumulative particle analysis by the HIAC showed the only formulation with significant sub-visible particle count was the “control” formulation, which showed an increasing trend with more aggressive stability conditions (higher temperatures for longer time). Other

formulations (M06-M08) showed only non-significant particle counts, independent of time or condition.

Functional Bioassay (Potency)

[0092] Functional FVIII bioassay results are presented below in **Table 16**. Results showed no significant changes in potency when samples are kept at either $\leq -65^{\circ}\text{C}$ or $2-8^{\circ}\text{C}$ conditions. At 25°C , there was a downward trend in potency for all formulations (M05-M08) after 4 weeks. At 40°C , there was a significant drop in potency for the Ca/Mg negative formulation after 3 days.

Table 16. Functional Bioassay (Relative Potency)

Formulation	Material	Timepoint					
		T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}
-70°C							
M05 - Control	M05	75	Not tested				81
M06 - Ca Negative	M06	90					81
M07 - Ca / Mg Negative	M07	67					68
M08 - 8.5% Sucrose	M08	74					74
5°C							
M05 - Control	M05	75	Not tested				73
M06 - Ca Negative	M06	90					87
M07 - Ca / Mg Negative	M07	67					71
M08 - 8.5% Sucrose	M08	74					76
25°C							
M05 - Control	M05	75	Not tested	53	Not tested	37	
M06 - Ca Negative	M06	90		66		63	
M07 - Ca / Mg Negative	M07	67		39		41	
M08 - 8.5% Sucrose	M08	74		55		45	
40°C							
M05 - Control	M05	75	51	Not tested			
M06 - Ca Negative	M06	90	60				
M07 - Ca / Mg Negative	M07	67	22				
M08 - 8.5% Sucrose	M08	74	56				

SEC

[0093] SEC results are presented below. Percent aggregate results are shown in **Table 17** below. No significant trend in aggregation was seen at $\leq -65^{\circ}\text{C}$, $2-8^{\circ}\text{C}$, or 25°C stability conditions. At 40°C , there was an increasing trend in aggregates for the Ca/Mg negative formulation, up to 7% HMWS after one week. The “calcium negative” formulation showed approximately 3% HMWS at the end of one week at 40°C . Other formulations did not show

significant growth in aggregation when held at $\leq -65^{\circ}\text{C}$ or $2-8^{\circ}\text{C}$ for up to 5 weeks, 25°C for up to four weeks, and 40°C for up to 1 week.

Table 17. Concentration of Aggregates by SEC (%HMWS)

Formulation	Material	Timepoint						
		T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}	T _{5W}
-70°C								
M05 – Control	M05	NMT 1.0	NT				NMT 1.0	NMT 1.0
M06 – Ca Negative	M06	NMT 1.0					NMT 1.0	NMT 1.0
M07 – Ca / Mg Negative	M07	NMT 1.0					NMT 1.0	NMT 1.0
M08 – 8.5% Sucrose	M08	NMT 1.0					NMT 1.0	NMT 1.0
5°C								
M05 – Control	M05	NMT 1.0	NT		NMT 1.0	NT	NMT 1.0	NMT 1.0
M06 – Ca Negative	M06	NMT 1.0			NMT 1.0		NMT 1.0	
M07 – Ca / Mg Negative	M07	NMT 1.0			NMT 1.0		NMT 1.0	
M08 – 8.5% Sucrose	M08	NMT 1.0			NMT 1.0		NMT 1.0	
25°C								
M05 – Control	M05	NMT 1.0	NT		NMT 1.0	NT	NMT 1.0	NT
M06 - Ca Negative	M06	NMT 1.0			NMT 1.0		NMT 1.0	
M07 - Ca / Mg Negative	M07	NMT 1.0			NMT 1.0		2.6	
M08 - 8.5% Sucrose	M08	NMT 1.0			NMT 1.0		NMT 1.0	
40°C								
M05 – Control	M05	NMT 1.0	NMT 1.0	NMT 1.0	NT			
M06 – Ca Negative	M06	NMT 1.0	3.1	2.6				
M07 – Ca / Mg Negative	M07	NMT 1.0	5.3	7.0				
M08 – 8.5% Sucrose	M08	NMT 1.0	1.9	NMT 1.0				

*NMT: no more than. NT: Not tested.

[0094] Viral particle results are shown below in **Table 18**. All formulations held at $\leq -65^{\circ}\text{C}$, $2-8^{\circ}\text{C}$, or 25°C conditions showed no significant change in viral particle titer. However, at 40°C , there was a significant drop in particle titer for the Ca/Mg negative formulation after 3 days and 1 week. Other formulations (M05-M06, M08) did not show a significant drop in viral particle titer when held at $\leq -65^{\circ}\text{C}$ or $2-8^{\circ}\text{C}$ for up to 5 weeks, 25°C for up to four weeks, and 40°C for up to 1 week.

Table 18. Viral Particle Titer by SEC (vp/mL)

Formulation	Material	Timepoint						
		T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}	T _{5W}
-70°C								
M05 – Control	M05	2.46E+13	NT				2.43E+13	2.54E+13
M06 – Ca Negative	M06	2.44E+13					2.41E+13	2.61E+13
M07 – Ca / Mg Negative	M07	2.02E+13					2.01E+13	2.11E+13
M08 – 8.5% Sucrose	M08	2.13E+13					2.16E+13	2.27E+13
5°C								
M05 – Control	M05	2.46E+13	NT		2.55E+13	NT	2.45E+13	2.51E+13
M06 – Ca Negative	M06	2.44E+13			2.54E+13		2.41E+13	2.62E+13
M07 – Ca / Mg Negative	M07	2.02E+13			2.11E+13		1.99E+13	2.14E+13
M08 – 8.5% Sucrose	M08	2.13E+13			2.25E+13		2.16E+13	2.34E+13
25°C								
M05 – Control	M05	2.46E+13	NT		2.49E+13	NT	2.35E+13	NT
M06 – Ca Negative	M06	2.44E+13			2.53E+13		2.39E+13	
M07 – Ca / Mg Negative	M07	2.02E+13			2.01E+13		1.90E+13	
M08 – 8.5% Sucrose	M08	2.13E+13			2.23E+13		2.11E+13	
40°C								
M05 – Control	M05	2.46E+13	2.23E+13	2.14E+13	NT			
M06 – Ca Negative	M06	2.44E+13	2.16E+13	2.09E+13				

Formulation	Material	Timepoint						
		T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}	T _{5W}
M07 – Ca / Mg Negative	M07	2.02E+13	1.62E+13	1.59E+13				
M08 – 8.5% Sucrose	M08	2.13E+13	2.09E+13	2.06E+13				

[0095] UV 260/280 ratio results are presented in **Table 19** below. For UV 260/280 ratio, there was a recorded downward trend for the Ca/Mg negative formulation when held at 40°C. UV 260/280 is considered a surrogate for empty/full particle ratio. Other formulations (M05, M06, and M08) did not show any significant change in UV 260/280 ratio at any condition or timepoint, as shown.

Table 19. UV 260/280 Ratio by SEC

Formulation	Material	Timepoint						
		T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}	T _{5W}
-70°C								
M05 - Control	M05	1.33	NT				1.33	1.33
M06 - Ca Negative	M06	1.32					1.32	1.33
M07 - Ca / Mg Negative	M07	1.32					1.32	1.33
M08 - 8.5% Sucrose	M08	1.32					1.32	1.33
5°C								
M05 - Control	M05	1.33	NT		1.33	NT	1.32	1.38
M06 - Ca Negative	M06	1.32			1.33		1.33	
M07 - Ca / Mg Negative	M07	1.32			1.33		1.32	1.33
M08 - 8.5% Sucrose	M08	1.32			1.33		1.32	1.33
25°C								
M05 - Control	M05	1.33	NT		1.34	NT	1.33	NT
M06 - Ca Negative	M06	1.32			1.32		1.33	
M07 - Ca / Mg Negative	M07	1.32			1.33		1.33	
M08 - 8.5% Sucrose	M08	1.32			1.32		1.32	
40°C								
M05 - Control	M05	1.33	1.33	1.33	NT			
M06 - Ca Negative	M06	1.32	1.32	1.32				
M07 - Ca / Mg Negative	M07	1.32	1.28	1.25				
M08 - 8.5% Sucrose	M08	1.32	1.33	1.33				

pH

[0096] There were no significant changes in pH at any condition and timepoint for all formulations, as shown in **Table 20** below.

Table 20. pH

Formulation	Material	Timepoint						
		T ₀	T _{3D}	T _{1W}	T _{2W}	T _{3W}	T _{4W}	T _{5W}
-70°C								
M05 - Control	M05	7.3	NT				7.2	7.3
M06 - Ca Negative	M06	7.3					7.3	7.3
M07 - Ca / Mg Negative	M07	7.3					7.3	7.4
M08 - 8.5% Sucrose	M08	7.2					7.3	7.3
5°C								
M05 - Control	M05	7.3	NT		7.3	NT	7.2	7.3
M06 - Ca Negative	M06	7.3			7.3		7.3	7.3
M07 - Ca / Mg Negative	M07	7.3			7.2		7.3	7.3
M08 - 8.5% Sucrose	M08	7.2			7.2		7.2	7.2
25°C								
M05 - Control	M05	7.3	NT		7.3	NT	7.2	NT
M06 - Ca Negative	M06	7.3			7.3		7.3	
M07 - Ca / Mg Negative	M07	7.3			7.3		7.2	
M08 - 8.5% Sucrose	M08	7.2			7.2		7.2	
40°C								
M05 - Control	M05	7.3	7.1	7.1	NT			
M06 - Ca Negative	M06	7.3	7.3	7.3				
M07 - Ca / Mg Negative	M07	7.3	7.3	7.3				
M08 - 8.5% Sucrose	M08	7.2	7.1	7.2				

Mass Spectrometry Characterization (Deamidation)

[0097] Deamidation results as determined by mass spectrometry are presented below in **Table 21**. Only samples stored at 40°C were used for this testing as samples stored at lower temperatures are not expected to show significant differences between formulations over the 5-week study period.

[0098] Results are presented using the T₀ of the control formulation as the T₀ of all formulations. Among all the deamidation sites, the N57G and N94H hotspots in AAV VP1 are the two which showed the most significant increase in deamidation over 3 weeks at 40°C. Other sites also showed an increase of deamidation, but to a much lesser extent (not shown). These two hotspots are especially important, as they are expected to affect transduction efficiency. The Ca/Mg negative formulation showed the most significant increase in deamidation. The other

formulations (M05, M06, and M08) all had similar levels of deamidation increase after 3 days and 3 weeks at 40°C.

Table 21. Deamidation by Mass Spectrometry (% at 40°C)

Formulation	Material	Timepoint					
		N57G (VP1)			N94H (VP1)		
		T ₀ *	T _{3D}	T _{3W}	T ₀ *	T _{3D}	T _{3W}
M05 - Control	M05	7.8	14.4	41.1	5.1	6.7	30.4
M06 - Ca Negative	M06	7.8	18.6	43.5	5.1	6.4	33.9
M07 - Ca / Mg Negative	M07	7.8	23.1	57.4	5.1	8.6	54.6
M08 - 8.5% Sucrose	M08	7.8	16.0	38.7	5.1	6.0	23.7

*T₀ was tested on M05 material and used as T₀ for all formulations.

Conclusions

[0099] Three possible variations on the current SB-525 formulation were evaluated: one with no calcium, one with neither calcium nor magnesium, and one with a higher concentration of sucrose. Examples 1 and 2 show the random generation of very many white-flake particles appearing in both buffer and drug product, especially after freeze/thaw cycling followed by being kept on stability at either 25°C for longer than 1 week, or 40°C for 3 days or longer. The particles were identified as calcium phosphate.

[0100] The results here demonstrate that calcium chloride is a non-essential excipient for product stability at intended and stressed storage conditions. Additionally, no significant difference was observed for “calcium negative” and “high sucrose” formulations for up to 2 weeks at 40°C, 4 weeks at 25°C, and 5 weeks at 5°C and -70°C. However, it was observed that removal of both calcium and magnesium from the formulation resulted in changes in drug product quality attributes, including showing larger drops in vg and vp titer as well as higher deamidation and aggregation, especially when kept at 40°C. In addition, it was found that through 5 weeks, only the previous (control) formulation continued to demonstrate the white flake-like particles.

Example 4: Further Reformulation of SB-525

[0101] In order to prevent particle generation, three types of changes to the formulation were tested in this Example: (1) increasing the sucrose concentration from 1% to 5-10% to increase stability against FT stress, (2) removing the divalent cations to remove the source of the particles

and prevent potential solubility issues, and (3) increase NaCl to shift the ionic strength of the formulation.

[0102] More specifically, the experiments in this Example assessed the stability of SB-525 buffer when reformulated to one of several variants in comparison to the previous (control) formulation; that is, 172 mM NaCl, 8.10 mM Na₂HPO₄, 2.68 mM KCl, 1.47 mM KH₂PO₄, 0.90 mM CaCl₂, 0.49 mM MgCl₂, 1% (w/v) sucrose, and 0.05% (w/v) poloxamer 188, pH 7.0-7.6. Formulations were filled at 5 mL into 6 mL AT vials, subjected to 5 uncontrolled FT cycles ($\leq$ -65°C to ambient), and then placed on thermal stability at 2-8°C (accelerated; liquid storage) and 25°C (stressed; liquid storage) to evaluate the various buffers against particle formation.

Study Design

[0103] Formulations tested, as shown in **Table 22**, were formulated by addition of sodium and potassium chlorides and phosphates to water, followed by addition of sucrose. If present, magnesium chloride, followed by calcium chloride, were then prepared in separate 50 mL solutions, transferred to the larger solution. Poloxamer 188 was then added to the solution. The formulations were then QS'd (quantum satis), pH tested, and filtered over a 0.22 μ m PES filter.

Table 22. Formulations

Material	Formulation	Difference from Control
0139-M03	200 mM NaCl	Increase NaCl to 200 mM
0139-M04	5% sucrose	Increase sucrose to 5%
0139-M05	8.5% sucrose	Increase sucrose to 8.5%
0139-M06	10% sucrose	Increase sucrose to 10%
0139-M07	No CaCl ₂	Removed CaCl ₂
0139-M08	No MgCl ₂ or CaCl ₂	Removed CaCl ₂ and MgCl ₂
0139-M09	Control	N/A

[0104] With these formulations, the AT vials were filled at 5 mL, subjected to 5 uncontrolled freeze/thaw cycles, and then placed on either 5°C or 25°C stability. The stability pull schedule is outlined in **Table 23**. In the Table, X refers to appearance (liquid); A refers to osmolality (freezing point depression), conductivity, viscosity, and density; B refers to light obscuration with HIAC; C refers to P188 concentration; and D refers to pH.

Table 23. Stability Schedule 0139-M03 to -M08

Temp (°C)	T ₀	FT1	FT2	FT3	FT4	FT5*	T _{5D}	T _{8D}	T _{2W}	T _{3W}	T _{4W}	T _{6W}	T _{8W}
5 (± 3)	XABD	X	X	XBD	X	XBD	X	X	XC	XB	XB	XB	XB
25 (± 3)							XBC		X				

* The 5 FT sample was the T₀ for the thermal stability samples.

Formulation Characterization

[0105] SB-525 buffer formulations were tested to determine conductivity, osmolality, density, and viscosity. P188 concentration was also measured at 2.5 weeks of 5°C storage, and 5 days of 25°C. The results to these tests are shown in **Table 24** below. The results in conductivity, osmolality, density, and viscosity for M03-M09 all follow expected trends – higher amounts of sucrose resulted in a higher osmolality, density, and viscosity; higher NaCl resulted in a higher osmolality and conductivity. Results for the concentration of P188 were within method variability for each set of results, confirming the correct P188 concentration was achieved during formulation.

Table 24. Formulation Characterization

Formulation	Conductivity (ms/cm)	Osmolality (mOsm/kg)	Density ¹ (g/mL)	Viscosity ² (Cp)	P188 (µg/mL)	
					T _{2.5W} 5°C	T _{5D} 25°C
M03 – 200 mM NaCl	20.20	433	1.011915	1.165	590	484
M04 – 5% Sucrose	16.12	515	1.026105	1.285	630	467
M05 – 8.5% Sucrose	15.37	647	1.039606	1.382	635	486
M06 – 10% Sucrose	15.12	704	1.045258	1.429	633	457
M07 – Ca Negative	17.91	380	1.010685	1.112	663	488
M08 – Ca / Mg Negative	17.88	379	1.010806	1.123	672	471
M09 – Control	17.86	384	1.010854	1.105	668	475

¹Temperature during density measurement was recorded as 20.0°C.

²Temperature during viscosity measurement was recorded as 20.0°C and shear rate was recorded as 2000 s⁻¹.

pH

[0106] No significant change in pH was observed for any formulation during FT cycling, as shown below in **Table 25**.

Table 25. pH Results

Formulation	FT Cycle		
	T ₀	FT3	FT5
M03 – 200 mM NaCl	7.2	7.3	7.2
M04 – 5% Sucrose	7.2	7.2	7.2
M05 – 8.5% Sucrose	7.2	7.2	7.2
M06 – 10% Sucrose	7.1	7.2	7.2
M07 – Ca Negative	7.3	7.3	7.3
M08 – Ca and Mg Negative	7.3	7.3	7.3
M09 – Control	7.2	7.3	7.2

Visible Appearance

[0107] All samples displayed B9 color and less than or equal to opalescence reference 1 for each condition and timepoint tested. Visible appearance results showed that only the control formulation (M09) consistently displayed too many to count (TMTC), white flake-like particles from 5 days to up to 8 weeks at 5°C and 25°C. The 5% sucrose formulation exhibited random generation of TMTC visible particles. Many samples also showed 1-5 fibrous particles. While these particles were not identified for this study, fiber-like particles are more characteristic of exogenous particles (e.g., filter particles) versus the flake-like particles that are characteristic of calcium phosphate. Also, appearance of fiber-like/exogenous particles was negligible for lab-based development studies where DP is not manufactured under tightly controlled environmental conditions.

Sub-Visible Particle Count

[0108] Sub-visible particle analysis by the HIAC showed the only formulation with significant and consistent sub-visible particle count was the control formulation (M09), exacerbated by ongoing storage at a higher stability condition (25°C). Other formulations (M03-M08) showed non-significant particle counts, independent of time or condition, and demonstrated no overall trend.

Conclusions

[0109] This Example evaluated the effects of freeze-thaw stress and thermal stability on several formulation variants from the previous (control) SB-525 formulation. Seven different formulations were filled into AT crystal closed vials, freeze-thaw cycled, and stored at 5°C and 25°C for 8 weeks. Across all study conditions, only the control formulation continued to demonstrate the white flake-like particles consistently.

[0110] After FT cycling, no formulation exhibited significant visible particles formation or changes in pH, and thus all formulations were placed on short-term stability at 5°C and 25°C over a period of 8 weeks. Upon placing on stability, the only formulation shown to consistently generate significant visible particles was the control formulation (M09). Although the particle results were comparable (i.e., no significant visible particle generation), the “200 mM NaCl” (M03), “5% sucrose” (M04), and “10% sucrose” (M06) formulations were not selected for continued reformulation studies. Particularly, the “5% sucrose” formulation (M04) exhibited some formation of white flake-like particles, indicating that this concentration of sucrose may not be sufficient to prevent particle generation due to FT stress, and “10% sucrose” (M06) would introduce manufacturing challenges due to its higher viscosity. Therefore, based on these results, the “8.5% sucrose” (M05), “calcium negative” (M07), and “calcium and magnesium negative” (M08) formulations appear promising.

Example 5: Final Formulation for SB-525

[0111] The composition and description of SB-525 (PF-07055480) Drug Substance (DS) is shown in **Table 26**. The formulation had a final pH of 7.3 ± 0.3 . Excipient concentration was obtained wholly or in part through the addition of the base buffer (Phosphate Buffered Saline). DS concentration target was $1.0\text{E}+13$ vg/mL ($0.5\text{E}+13 - 2.5\text{E}+13$ vg/mL), with a DS range of 50-250% of DS target.

Table 26. SB-525 Drug Substance Formulation Description

Raw Material	Chemical Formula	FW	Concentration (mM)	Concentration (g/L)
SB-525	N/A	N/A	$1.0\text{E}+13$ vg/mL	$1.0\text{E}+13$ vg/mL
Magnesium Chloride Hexahydrate	$\text{MgCl}_2 \bullet 6\text{H}_2\text{O}$	203.30	0.49	0.10
Potassium Chloride	KCl	74.5513	2.68	0.20

Raw Material	Chemical Formula	FW	Concentration (mM)	Concentration (g/L)
Potassium Phosphate Monobasic	KH ₂ PO ₄	136.086	1.47	0.20
Sucrose	C ₁₂ H ₂₂ O ₁₁	342.2965	1% (w/v)	10.00
Sodium Chloride	NaCl	58.44	172	10.05
Poloxamer 188 (P188)	(C ₃ H ₆ O • C ₂ H ₄ O) _x	8350	0.05% (w/v)	0.50
Sodium Phosphate Dibasic Anhydrous	Na ₂ HPO ₄	141.96	8.10	1.15
WFI	H ₂ O	18.015	q.s. to 1 mL	q.s. to 1 mL

[0112] The composition and description of SB-525 (PF-07055480) Drug Product is shown in **Table 27**. The DP had a final pH of 7.3 ± 0.3 . Excipient concentration was obtained wholly or in part through the addition of the base buffer (Phosphate Buffered Saline). DP concentration target was $1.0\text{E}+13$ vg/mL ($0.3\text{E}+13$ – $3.0\text{E}+13$ vg/mL), with a DP range of 30-300% of DP target.

Table 27. SB-525 Drug Product Formulation Description

Name of Ingredients	Function	Concentration (mM)	Unit Formula (mg/mL)	Unit Formula (mg/vial)
SB-525	Active Ingredient	$1.0\text{E}+13$ vg/mL	$1.0\text{E}+13$ vg/mL	$6.0\text{E}+13$ vg/vial
Sodium Phosphate Dibasic Anhydrous	Buffering Components	8.10	1.15	6.75
Potassium Phosphate Monobasic		1.47	0.20	1.20
Sodium Chloride	Stabilizers	172	10.05	60.30
Potassium Chloride		2.68	0.20	1.20
Magnesium Chloride Hexahydrate		0.49	0.10	0.60
Sucrose	Cryoprotectant	1% (w/v)	10.00	60.00
Kolliphor (Poloxamer) P 188	Surfactant	0.05% (w/v)	0.50	3.00
WFI	Solvent	q.s. to 1 mL	q.s. to 1 mL	q.s. to nominal fill volume

Container and Fill Volume

[0113] Filtered Drug Substance was packaged in sterile high-density polyethylene (HDPE) bottles. The DP was filled into AT 10 mL cyclo-olefin copolymer (COC) vials with in-place thermoplastic elastomer stoppers (Aseptic Technologies VIA-101800). Both DS and DP were

stored at -60°C to -90°C. Each vial contained 6.4 mL of DP (intended extractable volume 6 mL), with 6.4E+13 (nominal 6.0E+13) vg SB-525 AAV.

Characterization of DP

[0114] The viscosity, osmolality, density, and conductivity of the SB-525 DP were measured. The density was 1.0106 g/mL (20°C). The viscosity was 1.112 cP (20°C). The osmolality was 374 mOsm/kg. The conductivity was 17.80 mS/cm (20°C).

Example 6: AAV6 Viral Vector Drug Product Formulations

[0115] Examples 6-8 below describe an additional reformulation and accelerated stability study performed to evaluate AAV2/6 viral vector drug product formulations containing various concentrations of divalent cation salts. As in the previous example, the viral vector has an AAV6 capsid and a recombinant genome containing AAV2 ITRs. The viral vector here carries a transgene that encodes alpha-L-iduronidase (IDUA). This transgene would help patients who is IDUA-deficient, such as patients with mucopolysaccharidosis type I (MPS I), also known as Hurler Syndrome (a lysosomal storage disorder). The AAV genome is shown as sequence No. 28 in Table 5 of US2020/0246486. Because the viral genome is located within the AAV capsid and not exposed to the product formulation, the observations made in Examples 6-8, like those in the previous examples, are expected to be applicable to AAV6 vectors that carry other transgenes.

[0116] The study also was performed to identify a new product formulation that does not form precipitates following freeze/thaw cycles. The stability of four AAV6 viral vector drug product formulations containing various concentrations of divalent cation salts was assessed.

[0117] The formulations tested contained various concentrations of divalent cation salts (Ca^{2+} and Mg^{2+}). As described above, prior AAV6 formulations as well as formulation buffer without AAV6 were all observed to form small amounts of visible precipitates following multiple freeze/thaw cycles. The precipitated particles were determined to be comprised of calcium phosphate salts. Since the drug product will experience freeze/thaw cycles during typical use, the chance of salt precipitation would pose a risk to patient safety as well as a risk to product quality. The study shown in Examples 6-8 below aimed to find improved formulations that do not have such a precipitation issue.

[0118] In the study, samples of AAV6 viral vector product were prepared in different product formulations by tangential flow filtration and spiking samples with higher concentration “spike buffers.” The reformulated samples were then incubated at various accelerated stability conditions, and the samples’ quality attributes were evaluated following incubation. Quality attribute results were assessed with consideration for known method variability where appropriate. Outliers and trending results over the study duration were identified for each of the four formulations. A scoring of Pass, Neutral, or Fail was assigned for each quality attribute and formulation based on the results at each study endpoint. Finally, the overall quality attribute stability scorings of each four formulations were compared against each other. The study is described in detail below.

Sample Preparation

[0119] AAV6 product formulations prepared herein contained AAV6 donor vector at approximately $1.0\text{E}+13$ vg/mL suspended in formulation buffer F0 as described in Table 1 below. Buffer F0 is based upon a recipe of Dulbecco’s phosphate buffered saline that includes divalent cations calcium (as CaCl_2) and magnesium (as MgCl_2), approximately 35 mM additional sodium chloride, and is formulated to 1% w/v sucrose and 0.05% w/v poloxamer 188.

[0120] In order to transfer the AAV6 material into alternate buffer formulations, buffer exchange by tangential flow filtration was performed as shown in **FIG. 2**. Briefly, several vials of the AAV6 drug product formulated in buffer F0 were thawed and transferred to a single Ultrafiltration/Diafiltration (UF/DF) reservoir. Ultrafiltration was first performed to concentrate the material approximately 2x. Then the concentrated material was diafiltered with at least 10x equivalent volumes of formulation buffer F1. The recovered intermediate material, labeled “A” in **FIG. 2**, should thus contain approximately $2.0\text{E}+13$ vg copies/mL suspended in formulation buffer F1.

[0121] The intermediate UF/DF product “A” was then separated into three volumes and formulated into three different buffers. Three additional formulation buffers (F1, F2, and F3; **Table 28**) were derived from Dulbecco’s phosphate buffered saline with additional sodium chloride, sucrose, and poloxamer 188 (e.g., Pluronic® F-68). These three buffers contained no calcium and instead had different amounts of magnesium. Buffer F1 contained no calcium or magnesium components, Buffer F2 contained the equivalent molar concentration of magnesium

as F0, and Buffer F3 contained additional magnesium accounting for the moles of calcium removed from buffer F0.

Table 28. Composition of Formulation Buffers

Component (unit)	F0	F1	F2 (6x conc)	F2 (1x conc)	F3 (6x conc)	F3 (1x conc)
NaCl (mM)	171.8	171.8	171.8	171.8	171.8	171.8
KCl (mM)	2.7	2.7	2.7	2.7	2.7	2.7
Na ₂ HPO ₄ (mM)	8.1	8.1	8.1	8.1	8.1	8.1
KH ₂ PO ₄ (mM)	1.5	1.5	1.5	1.5	1.5	1.5
CaCl ₂ (mM)	0.9	0.0	0.0	0.0	0.0	0.0
MgCl ₂ (mM)	0.5	0.0	3.0	0.5	8.4	1.4
Sucrose (% w/v)	1.00%	1.00%	1.00%	1.00%	1.00%	1.00%
Poloxamer 188 (% w/v)	0.05%	0.05%	0.05%	0.05%	0.05%	0.05%
Target pH	7.1-7.5					

[0122] To simplify sample preparation and to allow preparation of different buffers, 6-times higher concentrated “spike buffers” of buffers F2 and F3 were prepared. As shown in **FIG. 3**, these F2 and F3 “6x conc” buffers when mixed at a 1:5 ratio with F1 yielded the intended composition of buffers F2 and F3, respectively.

[0123] To prepare these formulation buffers, individual dry components were weighed off using a microbalance. The components were mixed and dissolved with deionized water for injection. These buffers were then titrated with 0.1 M NaOH or 0.1 M HCl to the intended pH range if needed. Finally, these buffers were formulated with sucrose and poloxamer 188 to the desired concentrations.

[0124] The compositions of the prepared drug product formulations were measured with several semiquantitative assays as a check, the results of which are reported in **Table 29**. Off-the-shelf colorimetric assays were used to measure concentrations of Ca²⁺ (Bio Vision Cat# K380-250) and Mg²⁺ (Bio Vision Cat#K385-1 00). For both colorimetric assays, the sample concentrations were within the assays’ dynamic linear ranges, and the reported signals and calculated concentrations trended well with the theoretical sample compositions. Critically, the calcium signals for formulations F1, F2, and F3 all matched the assay negative controls. Thus, these samples contained no calcium, and the TFF sample preparation step was effective at removing any calcium that was present in the F0 starting material. The magnesium assay also

confirmed Mg^{2+} was also removed by TFF and spiked into the appropriate concentrations in F2 and F3.

Table 29. Experimentally Measured Parameters of Formulated Drug Products

Component (unit)	F0	F1	F2 (1x conc)	F3 (1x conc)
CaCl_2 (mM, by colorimetric assay)	1.02	0.07	0.09	0.1
MgCl_2 (mM, by colorimetric assay)	0.29	-0.04	0.38	0.9
VG titer (vg/ml, by qPCR)	1.54E+13	1.77E+13	1.91E+13	1.87 E+13
Particle Concentration (p/mL, by MADLS)	7.79E+12	6.38E+12	7.30E+12	8.23E+12

*MADLS: Multi-Angle Dynamic Light Scattering.

[0125] qPCR was also performed to measure vg titer, and multi-angle dynamic light scattering (MADLS) was performed to measure particle concentration. Again, these results were in the expected ranges from sample preparation, indicating that the buffer exchange and formulation steps performed as intended.

[0126] The bulk formulated products were sterile filtered and filled at 0.5 mL/vial into 2 mL cyclic olefin polymer (COP) vials (West Pharma/Daikyo CZ Cat # 19550057).

Sample Incubation

[0127] Vial drug products were incubated at various conditions described below.

Freeze/thaw cycling was performed with pull dates listed in **Table 30**.

Table 30. Freeze/Thaw Cycling Sample Pull Dates

Condition	Start Date	3X Cycles	5X Cycles	10X Cycles
Freeze/Thaw cycling -80°C/Ambient	D0	D2	D4	D11

[0128] Samples were subject to “accelerated, stressed” conditions (40°C/75% relative humidity (RH)) or to “accelerated, ambient” conditions (25°C/60% RH), with pull dates at D3, D7, 2 weeks, 1 month, 2 months, 3 months, or 6 months. Samples were stored at 2-8°C and analyzed within 24 hours of their pull time.

[0129] All pulled samples were aliquoted into polycarbonate flip-top centrifuge tubes for the intended sample testing. Aliquots were stored at 2-8°C during active testing for no longer than 2 months. If testing was not performed within a reasonable time from the pull date, the aliquots

were frozen at $\leq -65^{\circ}\text{C}$ until testing could occur. All retain samples and remaining materials following testing were frozen at $\leq -65^{\circ}\text{C}$ for long term storage.

Example 7: Testing of AAV6 Drug Product Formulations

[0130] Sample testing performed included the quality attributes listed in **Table 31**.

Abbreviations are shown as below: DLS, dynamic light scattering; MADLS, multi-angle dynamic light scattering; HMWS, high molecule weight species; SE-LC, size exclusion liquid chromatography; and AEX-LC, anion exchange liquid chromatography. Appearance of the solutions was assessed to confirm that they were clear, colorless, and free of particulates. A solution that was hazy, cloudy, or containing visible particles did not meet the criteria.

Table 31. Product Quality Attributes Evaluated

Category	Quality Attribute	Testing Methods
General	Appearance, Liquid	Visual inspection
	pH	Potentiometric
Purity	Monomer Peak Mean Size	DLS
	Polydispersity Index	DLS
	Monomer Fraction	DLS
	HMWS Fraction	DLS
	HMWS Fraction	SE-LC
	% Full Capsids	AEX-LC
Strength	Vector Genome Titer	qPCR
	Infectious Titer	TCID ₅₀
	VG/TCID ₅₀ Ratio	Calculation
	Capsid Concentration	ELISA
	Monomer Concentration	SE-LC
	Monomer Concentration	MADLS

Initial Samples

[0131] Quality attribute results from all initial formulated drug product samples are shown in **Table 32** and **Table 33**. Results are grouped by the quality attribute categories of general, purity, and strength.

Table 32. General and Purity Quality Attribute Results of Initial (to) Samples

Formulation	F0	F1	F2	F3
Appearance (Pass/Fail)	Pass	Pass	Pass	Pass
pH	7.3	7.3	7.2	7.2

Monomer Peak Mean Size (nm)	29.7	33.0	33.6	33.4
Polydispersity Index	0.07	0.12	0.14	0.12
% Monomer, DLS	100.0%	100.0%	100.0%	100.0%
% HMWS, DLS	0.00%	0.00%	0.00%	0.00%
% HMWS, SE-LC	0.2%	0.5%	0.5%	0.4%
% Full Capsids	89.2%	89.1%	89.3%	89.1%

Table 33. Strength Quality Attribute Results of Initial (t₀) Samples

Formulation	F0	F1	F2	F3
VG Titer (vg/mL)	7.9E+12	8.3E+12	9.4E+12	7.7E+12
Infectious Titer (TCID ₅₀ /mL)	9.3E+09	9.3E+09	1.1E+10	1.1E+09
VG/TCID ₅₀	1200	1100	1100	8400
Capsid Titer (capsids/mL)	5.0E+12	5.5E+12	5.2E+12	6.7E+12
Monomer Conc, SE-LC (particles/mL)	9.9E+12	9.5E+12	9.9E+12	9.7E+12
Monomer Conc, MADLS (particles/mL)	7.8E+12	6.4E+12	7.3E+12	8.2E+12

[0132] The initial infectious titer result was usually low for F3 when compared against the other initial samples. This also resulted in a much higher vg/TCID ratio than the other samples. This was attributed primarily to the variability of the TCID₅₀ assay. None of the other results appeared to be substantially affected.

Freeze/Thaw Cycled Samples

[0133] Test results from freeze/thaw cycling of the four different formulations are shown in **FIGs. 4A-D**. The data show that the freeze/thaw cycling did not affect solution pH, monomer capsid size, sample size distribution, or full capsid content; nor did it generate significant amounts of HMWS or soluble aggregates in all four formulations. Freeze/thaw cycling did, however, lead to visible appearance failures in the initial F0 control formulations. Following 5x F/T cycles, small white flecks of precipitates could be observed, and these were characteristic of the calcium phosphate precipitates previously observed. These precipitates were not detected in formulations F1, F2, or F3. The data also show that freeze/thaw cycling also did not appear to significantly affect any product strength attributes when the variability of test methods is also considered.

25°C/60% RH (Ambient) Accelerated Samples

[0134] Test results from ambient condition (25°C/60% RH) incubation of the four different formulations are shown in **FIGs. 5A-D**. The data show that appearance and pH were not

affected by the 25°C condition. The sample size distribution appeared to broaden, increasing the polydispersity index (PDI), and the monomer peak mean size grew slightly as the incubation time increased. High weight molecular species (HMWS) also began to form over incubation time. In general, these HMWS were detected more readily by SE-LC than by DLS. Charge separation of full and empty capsids appears to have been lost following 1-2 months of incubation at 25°C, thus the % full capsids could not be measured by AEX-LC. Capsid titer and monomer concentration were not affected by the 25°C condition, nor did the vg titer appear to be affected. TCID₅₀ decreased at the 3-month and 6-month time point for F0, F1, and F2 but the same result was not observed with F3.

40°C/75% RH Stressed Samples

[0135] Test results from stressed condition (40°C/75% RH) incubation of the four different formulations are shown in **FIGs. 6A-D**. The data show that the 40°C/75% RH condition had a much greater impact on several product quality attributes. The particle size distribution broadened substantially and HMWS were forming. SE-LC appeared to be more sensitive to detecting smaller amounts of HMWS than DLS, but when large amounts of HMWS (and possibly larger molecular weight species) were present, these particles might not have been detected by SE-LC. It is likely that large molecular weight species were being filtered or not passing through the LC column to reach the detector in these cases.

[0136] Following 3 months at 40°C, all solutions began to appear hazy, possibly due to the formation of large insoluble aggregates or HMWS. The solution pH did not appear to be affected. TCID₅₀ did decrease substantially in all samples after 1 month at 40°C/75% RH with the smallest reduction from initial seen in F3.

Example 8: Study Endpoint Assessment

[0137] An intra-sample analysis was performed for each formulation. Within each formulation, each study endpoint result for each quality attribute was scored against the acceptance and failure criteria as described in **FIG. 7**. The intended stability study endpoints were as follows: 10x freeze/thaw cycles of no less than 12 hours at $\leq -65^{\circ}\text{C}$ and no less than 6 hours at ambient temperature per cycle, 3 months at 25°C/60% RH, and 1 month at 40°C/75% RH. At each study endpoint, the quality attributes as listed in **FIG. 7** were determined based on

known and established method variability as well as a holistic assessment of all data generated in this study.

[0138] To perform the inter-sample analysis or a relative assessment of each formulation, each endpoint score for each quality attribute was compared by formulation. When considering a single quality attribute, an overall Pass, Neutral, or Fail scoring was assigned for each formulation by holistically comparing the relative endpoint scorings across the four considered formulations. For example, if all endpoint results had passed acceptance criteria for a given formulation, an overall Pass score was assigned for that quality attribute for that formulation. The number of neutral and failing endpoint conditions were considered when assigning an overall quality attribute scoring for that formulation. For certain quality attributes such as appearance, one failing endpoint score was enough to score an overall Fail for that formulation. For other quality attributes such as monomer concentration, one failing endpoint at the 40°C condition was encountered for all formulations (except F3); thus, this failing endpoint was not scored as severely. TCID₅₀ results were included for information only and a scoring was not assigned.

[0139] These intra-sample analysis scores are shown in **FIG. 8** and **FIG. 9**. A summary table of all overall quality attribute scores for each formulation is shown in **Table 34**, where formulation F0 is used as a control.

Table 34. Inter-Sample Formulation Assessment Summary

Test (method)	F0	F1	F2	F3
Appearance, Liquid (Visual Inspect)	Fail	Pass	Pass	Pass
pH (Potentiometric)	Pass	Pass	Pass	Pass
Size Distribution (DLS)	Neutral	Fail	Neutral	Neutral
HMWS Fraction (DLS)	Pass	Fail	Pass	Pass
HMWS Fraction (SE-LC)	Neutral	Fail	Neutral	Pass
% Full Capsids (AEX-LC)	Pass	Pass	Pass	Pass
VG Titer (qPCR)	Pass	Fail	Pass	Neutral
Capsid Conc. (ELISA)	Pass	Neutral	Pass	Pass
Monomer Conc (SE-LC)	Neutral	Neutral	Neutral	Pass
Monomer Conc. (MADLS)	Neutral	Neutral	Neutral	Pass

[0140] The overall quality attribute scorings over the duration of the accelerated stability study demonstrate that Formulations F3 and F2 are superior to F0 and F1. Formulation F3 appears to have performed the best, with eight overall passing quality attribute scores and only two neutral scores. Formulation F2 also performed well. F0 had an overall failing score for appearance, confirming the initial goal for improving this product formulation. F1 appeared to have the greatest amount of aggregation with the greatest amount of HMWS detected.

Example 9: Extended Temperature Tests

[0141] Stability of SB-525 (PF-07055480) drug product provided at 1.00E+13 vg/mL in formulation buffer was tested at extended temperatures. The formulation buffer contained the following ingredients: 0.49 mM MgCl₂, 2.68 mM KCl, 1.47 mM KH₂PO₄, 172 mM NaCl, 8.10 mM Na₂HPO₄, 1% (w/v) sucrose, 0.05% (w/v) poloxamer 188, pH 7.4.

[0142] Vials containing the drug product were stored at (i) -70°C for 0 days (T₀), (ii) at -150°C for 3 days (T_{3-Days}), and (iii) at -150°C for 14 days (T_{14-Days}). Drug product samples were then analyzed for (i) appearance (liquid), (ii) reduced CE-SDS, (iii) SEC titer_{260/280}, and (iv) *in vitro* FVIII activity (bioassay). The formulation buffer sample was also analyzed for (i) container closure integrity.

[0143] As shown in **Table 35** below, there was no significant change in color, clarity, reduced CE-SDS, SEC titer_{260/280} (a surrogate measure for empty:full capsid ratio), or biological activity (% relative potency) when the samples were stored at the tested conditions including storage at -150°C for up to fourteen days. No impact on container closure integrity was observed after fourteen days storage at -150°C.

Table 35. Stability Results

Test		T ₀	T _{3-Days}	T _{14-Days}
Appearance, Liquid	Color	<=B9	<=B9	<=B9
	Clarity	<= Ref. I	<= Ref. I	<= Ref. I
	Visible Particles	EFVP	EFVP	EFVP
Container Closure Integrity		N.T.	N.T.	Pass
Capsid Protein Purity (CE): % Capsid		95.0	94.5	94.9

Test	T ₀	T _{3-Days}	T _{14-Days}
Capsid Protein Purity (CE): % Other	5.0	5.5	5.1
Size Distribution (SEC): % Monomer	98.8	99.0	99.5
Size Distribution (SEC): % Aggregates	1.2	NMT 1.0	NMT 1.0
Particle Content (SEC): 260/280 Ratio	1.34	1.34	1.33
<i>in vitro</i> FVIII Activity/ Functional Bioassay (% Relative Potency)	67%	69%	70%

EFVP: Essentially Free of Visible Particles

TMTC: Too many to Count

N.T.: Not tested

Example 10: Shipping Simulation and Extended Temperature Tests

[0144] Stability of the SB-525 drug product under simulated shipping conditions, such as shock, pressure, drop, and vibration, was tested. The SB-525 drug product was provided at 1.00E+13 vg/mL in formulation buffer containing the following ingredients: 0.49 mM MgCl₂, 2.68 mM KCl, 1.47 mM KH₂PO₄, 172 mM NaCl, 8.10 mM Na₂HPO₄, 1% (w/v) sucrose, 0.05% (w/v) poloxamer 188, pH 7.4.

[0145] Test sample vials containing the drug product underwent concurrently applied transport hazards (e.g., shock, pressure, drop, and vibration), utilizing a worst-case global transportation profile. Control sample vials were not subjected to those concurrently applied transport hazards. As part of the simulation, test vials were held for 40 hours at -35°C and subsequently for another 40 hours at -70°C. After simulation, control and test drug product vials were stored at -70°C and tested at 0 days (T₀), 6 months (T_{6-Months}), and 10 months (T_{10-Months}). Formulation buffer vials were stored at -70°C and tested at T₀. Drug product samples were analyzed for (i) appearance (liquid), (ii) reduced CE-SDS, (iii) SEC titer_{260/280}, and (iv) *in vitro* FVIII activity (bioassay). Formulation buffer samples were analyzed for (i) container closure integrity.

[0146] As shown in **Table 36** below, there was no significant change in color, clarity, reduced CE-SDS, SEC titer_{260/280} (a surrogate measure for empty:full capsid ratio), or biological activity (% relative potency) when the samples were stored at the tested conditions for up to ten months. There was no impact on container closure integrity at T₀.

Table 36. Stability Results

Test		Sample	T ₀	T _{6-Months}	T _{10-Months}
Appearance, Liquid	Color	Control	<=B9	<=B9	<=B9
		Test	<=B9	<=B9	<=B9
	Clarity	Control	<= Ref. I	<= Ref. I	<= Ref. I
		Test	<= Ref. I	<= Ref. I	<= Ref. I
	Visible Particles	Control	EFVP	EFVP	EFVP
		Test	EFVP	EFVP	EFVP
Container Closure Integrity		Control	Pass	N.T.	N.T.
		Test	Pass	N.T.	N.T.
Genome Titer (qPCR)		Control	1.11E+13	1.02E+13	7.96E+12
		Test	9.07E+12	9.29E+12	7.81E+12
Capsid Titer (ELISA)		Control	1.71E+13	1.28E+13	1.30E+13
		Test	1.69E+13	1.19E+13	1.34E+13
Capsid Protein Purity (CE): % Capsid		Control	94.0	94.3	95.0
		Test	94.5	94.6	94.7
Capsid Protein Purity (CE): % Other		Control	6.0	5.7	5.0
		Test	5.5	5.4	5.3
Size Distribution (SEC): % Monomer		Control	99.2	99.1	99.3
		Test	99.3	98.8	99.4
Size Distribution (SEC): % Aggregates		Control	NMT 1.0	NMT 1.0	NMT 1.0
		Test	NMT 1.0	1.2	NMT 1.0
Particle Content (SEC): 260/280 Ratio		Control	1.33	1.32	1.34
		Test	1.33	1.32	1.33
<i>in vitro</i> FVIII Activity/ Functional Bioassay (% Relative Potency)		Control	72	58	62
		Test	63	58	62

Example 11: Drug Product Stability Data – 24 Months

[0147] The purpose of this study was to establish long term (24 months) stability of the SB-525 drug product at the intended storage temperature (-70°C). For this study, the SB-525 drug product (DP) was purified from Sf9 insect cells and formulated at a target of 1.00E+13 vector vg/mL in 8.10 mM Na₂HPO₄, 1.47 mM KH₂PO₄, 0.49 mM MgCl₂, 2.68 mM KCl, 172 mM NaCl, 1% (w/v) sucrose, and 0.05% (w/v) poloxamer P188, pH 7.3 ± 0.3. The DP was then filled at 6.4 mL into 10 mL Aseptic Technologies (AT) crystal closed vials, and stored at -60°C to -90°C.

[0148] Results from the DP stability study, including 5 cycles of uncontrolled F/T and 24-hour agitation, are shown in **FIG. 10**. The data show that there was no apparent trend observed in color, clarity, pH, subvisible cumulative particulate matter, capsid purity and VP ratio by R-CGE, %HMMS and UV 260/280 by SEC-HPLC, and capsid titer up to 24 months at -70°C, 12

months at 5°C, and 3 months at 25°C / 60% RH. The 12 month, 5°C sample reported “1 long fiber particle” in solution which appeared to be extrinsic in nature based on morphology and testing in uncontrolled environmental conditions (i.e., open lab space); all other samples reported “essentially free of visible particles.” No change in these quality attributes were observed after five cycles of uncontrolled F/T or 24-hour agitation.

[0149] For all the time points and conditions tested, the variability in the results for genome titer and infectious virus titer as well as infectivity ratio were within expected assay variability and did not demonstrate an apparent trend. The results for *in vitro* relative potency, a sensitive and precise stability-indicating method, show a slight downward trend in potency at the intended storage condition (-70°C) after 24 months, but the results are within the clinical stability acceptance criterion. A stability-indicating downward trend in potency at 3 months at 25°C / 60% RH and at 12 months at 5°C was also observed.

[0150] Additionally, Container Closure Integrity (CCIT) was tested at 18 months at the intended storage condition (-70°C) by a validated headspace analyzer and recorded a passing result. Six contingency vials were used for this analysis. P188 concentration was also stable for up to 18 months at -70°C, 12 months at 5°C, and 3 months at 25°C / 60% RH. A slight decrease in P188 concentration was observed at 24 months at -70°C.

CLAIMS

1. A pharmaceutical composition comprising
a recombinant adeno-associated virus (rAAV) vector,
sodium chloride (NaCl)
potassium chloride (KCl),
disodium phosphate (Na_2HPO_4),
monopotassium phosphate (KH_2PO_4),
magnesium chloride (MgCl_2),
a polyol, and
a poloxamer,
optionally wherein the composition comprises no more than about 0.1 mM calcium chloride and
has a pH of about 7.1 to about 7.5.
2. The composition of claim 1, wherein the polyol is sucrose.
3. The composition of claim 1 or 2, wherein the poloxamer is poloxamer 188.
4. The composition of any one of claims 1-3, wherein the composition contains about 0.1 to
about 2.0 mM magnesium chloride.
5. The composition of claim 4, wherein the composition contains about 0.5 mM or higher
magnesium chloride.
6. The composition of claim 4, wherein the composition contains about 1.4 mM magnesium
chloride.
7. The composition of claim 4, wherein the composition contains more than about 1.3 mM
magnesium chloride.

8. The composition of any one of claims 1-7, wherein the composition contains about 150 to about 200 mM, optionally about 172 mM, sodium chloride.
9. The composition of any one of claims 1-8, wherein the composition contains about 2.5 to about 3.0 mM, optionally about 2.7 mM, potassium chloride.
10. The composition of any one of claims 1-9, wherein the composition contains about 5 to about 10 mM, optionally about 8 mM, disodium phosphate.
11. The composition of any one of claims 1-10, wherein the composition contains about 1.0 to about 2.0 mM, optionally about 1.5 mM, monopotassium phosphate.
12. The composition of any one of claims 1-11, wherein the composition contains about 0.5% to about 2% (w/v), optionally about 1% (w/v), sucrose.
13. The composition of any one of claims 1-12, wherein the composition contains about 0.01% to about 0.1% (w/v), optionally about 0.05% (w/v), poloxamer 188.
14. A pharmaceutical composition comprising
 - a recombinant adeno-associated virus (rAAV) vector,
 - about 171.81 mM sodium chloride,
 - about 2.68 mM potassium chloride,
 - about 8.10 mM disodium phosphate,
 - about 1.47 mM monopotassium phosphate,
 - about 1.40 mM magnesium chloride,
 - about 1.00% (w/v) sucrose, and
 - about 0.05% (w/v) poloxamer 188,optionally wherein the composition comprises no more than about 0.1 mM calcium chloride and has a pH of about 7.1 to about 7.5.
15. A pharmaceutical composition comprising

a recombinant adeno-associated virus (rAAV) vector,
about 172 mM sodium chloride,
about 2.68 mM potassium chloride,
about 8.10 mM disodium phosphate,
about 1.47 mM monopotassium phosphate,
about 0.49 mM magnesium chloride,
about 1.00% (w/v) sucrose, and
about 0.05% (w/v) poloxamer 188,

optionally wherein the composition comprises no more than about 0.1 mM calcium chloride and has a pH of about 7.1 to about 7.5.

16. The composition of any one of claims 1-15, wherein the rAAV comprises a genome comprising an expression cassette for a therapeutic protein.

17. The composition of claim 16, wherein the therapeutic protein is a human factor VIII polypeptide.

18. A pharmaceutical composition comprising:

a recombinant adeno-associated virus (rAAV) vector, wherein the rAAV comprises a genome comprising an expression cassette for expressing a human factor VIII polypeptide,
about 171.81 mM sodium chloride,
about 2.68 mM potassium chloride,
about 8.10 mM disodium phosphate,
about 1.47 mM monopotassium phosphate,
about 1.40 mM magnesium chloride,
about 1.00% (w/v) sucrose, and
about 0.05% (w/v) poloxamer 188,

optionally wherein the composition comprises no more than about 0.1 mM calcium chloride and has a pH of about 7.1 to about 7.5.

19. A pharmaceutical composition comprising:

a recombinant adeno-associated virus (rAAV) vector, wherein the rAAV comprises a genome comprising an expression cassette for expressing a human factor VIII polypeptide

about 172 mM sodium chloride,
about 2.68 mM potassium chloride,
about 8.10 mM disodium phosphate,
about 1.47 mM monopotassium phosphate,
about 0.49 mM magnesium chloride,
about 1.00% (w/v) sucrose, and
about 0.05% (w/v) poloxamer 188,

optionally wherein the composition comprises no more than about 0.1 mM calcium chloride and has a pH of about 7.1 to about 7.5.

20. The composition of any one of claims 17-19, wherein the human factor VIII polypeptide comprises SEQ ID NO:1.

21. The composition of any one of claims 17-19, wherein the genome of the rAAV comprises SEQ ID NO:2 or nucleotides 131-5,024 of SEQ ID NO:2.

22. The composition of any one of claims 1-21, the composition contains the rAAV at about $1.0\text{E}+12$ to about $1.0\text{E}+14$ vector genomes (vg) per mL, optionally about $1.0\text{E}+13$ to about $5.0\text{E}+13$ vg per mL.

23. The composition of claim 22, wherein the composition contains about $1.0\text{E}+13$ vg per mL.

24. The composition of any one of the preceding claims, wherein the rAAV comprises an AAV6 capsid protein, optionally inverted terminal repeats (ITRs) of AAV2.

25. A vial comprising 5-10 mL, optionally 6.4 mL, of the composition of any one of claims 1-24.

26. The vial of claim 25, wherein the vial is made of cyclo-olefin copolymer.
27. The vial of claim 25 or 26, wherein the vial has an in-place thermoplastic elastomer stopper.
28. A method of treating a patient in need of a therapeutic protein, comprising administering to the patient the composition of any one of claims 1-24.
29. A method of increasing the serum level of factor VIII in a human subject in need thereof, comprising administering intravenously to the human subject the pharmaceutical composition of any one of claims 17-24, or the entire content of the vial of any one of claims 25-27.
30. The method of claim 29, wherein the human subject has hemophilia A.
31. A method of treating hemophilia A in a human subject in need thereof, comprising administering intravenously to the human subject the pharmaceutical composition of any one of claims 17-24, or the entire content of the vial of any one of claims 25-27.
32. Use of the pharmaceutical composition of any one of claims 1-24 for the manufacture of a medicament for treating a human subject in a method of any one of claims 28-31.
33. The pharmaceutical composition of any one of claims 1-24 for use in treating a human subject in a method of any one of claims 28-31.




Analysis	Sample	Result	Description
Digital imaging	F0		White particulate matter can be observed
Digital Imaging	F0 5x F/T		Large amounts of white particle matter can be observed
Digital Imaging	F0 after 1M at 25 °C		Large amounts of white particle matter can be observed

FIG. 1A


Polarized Light Microscopy	F0 after 1M at 25 °C	 Plane polarized  Crossed Polars	No birefringence seen under crossed-polars suggesting the material is either isotropic or amorphous
----------------------------------	----------------------	---	---

FIG. 1B



SEM/EDS	F0 5x F/T	<table><thead><tr><th></th><th>O</th><th>P</th><th>Ca</th></tr></thead><tbody><tr><td>Background square</td><td>28.94</td><td>32.76</td><td>0.00</td></tr><tr><td>Sample square 1</td><td>41.58</td><td>22.65</td><td>27.87</td></tr><tr><td>Sample square 2</td><td>41.55</td><td>21.14</td><td>24.13</td></tr></tbody></table>		O	P	Ca	Background square	28.94	32.76	0.00	Sample square 1	41.58	22.65	27.87	Sample square 2	41.55	21.14	24.13	There are large proportions of O, Ca and P in the sample
	O	P	Ca																
Background square	28.94	32.76	0.00																
Sample square 1	41.58	22.65	27.87																
Sample square 2	41.55	21.14	24.13																

FIG. 1C




FTIR	F0	 The FTIR spectrum of the F0 sample shows a broad absorption band centered around 3400 cm⁻¹, characteristic of O-H stretching in a salt. There are several sharp peaks in the fingerprint region, notably around 1650, 1550, 1450, 1350, 1250, 1150, and 1050 cm⁻¹. The x-axis is labeled 'Wavenumber (cm⁻¹)' and ranges from 4000 to 400. The y-axis is labeled 'Transmittance (%)' and ranges from 0 to 100.	FTIR suggest material is a salt from database search. F0, F0 5x F/T, and F0 1 M at 25°C samples all show similar spectra.
FTIR	F0 5x F/T	 The FTIR spectrum of the F0 5x F/T sample is very similar to the F0 sample, showing a broad O-H stretching band around 3400 cm⁻¹ and characteristic peaks in the fingerprint region at approximately 1650, 1550, 1450, 1350, 1250, 1150, and 1050 cm⁻¹. The axes are the same as the F0 spectrum.	FTIR suggest material is a salt from database search. F0, F0 5x F/T, and F0 1 M at 25°C samples all show similar spectra.
FTIR	F0 after 1M at 25°C	 The FTIR spectrum of the F0 sample after treatment with 1M at 25°C shows a broad O-H stretching band around 3400 cm⁻¹ and characteristic peaks in the fingerprint region at approximately 1650, 1550, 1450, 1350, 1250, 1150, and 1050 cm⁻¹. The axes are the same as the previous spectra.	FTIR suggest material is a salt from database search. F0, F0 5x F/T, and F0 1 M at 25°C samples all show similar spectra.

FIG. 1D




FIG. 1E



Digital Imaging	F0 5x F/T (lot 2)		White particle matter observed
FTIR	F0 5x F/T (lot 2)		Resembles spectra for F0 from previous lot
Raman microscopy	F0 5x F/T		Material identified as calcium phosphate

FIG. 1F




FIG. 2




FIG. 3




FIG. 4A




FIG. 4B

Formulation	Condition	Appearance	pH	Monomer Peak Mean Size (nm)	Polydisp. Index	% Monomer DLS	% HMWS, DLS	% HMWS, SE-LC	% Full Capsids
F0	Initial	Pass	7.3	29.7	0.07	100.0%	0%	0.2%	89.2%
	3x F/T	Pass	7.3	30.7	0.12	99.6%	0.4%	ND	ND
	5x F/T	Fail	7.3	30.4	0.08	100.0%	0%	ND	ND
	10x F/T	Fail	7.3	30.2	0.07	100.0%	0%	0.0%	89.3%
F1	Initial	Pass	7.3	33.0	0.12	100.0%	0%	0.5%	89.1%
	3x F/T	Pass	7.2	32.8	0.13	100.0%	0%	ND	ND
	5x F/T	Pass	7.2	34.3	0.13	100.0%	0%	ND	ND
	10x F/T	Pass	7.3	33.4	0.12	100.0%	0%	0.6 %	89.2%
F2	Initial	Pass	7.2	33.6	0.14	100.0%	0%	0.5%	89.3%
	3x F/T	Pass	7.2	34.0	0.14	100.0%	0%	ND	ND
	5x F/T	Pass	7.3	33.9	0.13	100.0%	0%	ND	ND
	10x F/T	Pass	7.2	33.6	0.13	100.0%	0%	0.5%	89.3%
F3	Initial	Pass	7.2	33.4	0.12	100.0%	0%	0.4%	89.1%
	3x F/T	Pass	7.2	33.3	0.12	100.0%	0%	ND	ND
	5x F/T	Pass	7.2	33.6	0.13	100.0%	0%	ND	ND
	10x F/T	Pass	7.3	33.7	0.13	100.0%	0%	0.3%	89.0%

ND: No data.

FIG. 4C

Formulation	Condition	VG Titer (vg/mL)	Infectious Titer (TCID ₅₀ /mL)	VG/TCID ₅₀	Capsid Titer (capsids/mL)	Monomer Conc, SE-LC (p/mL)	Monomer Conc, MADLS (p/mL)
F0	Initial	7.9E+12	9.3E+09	1200	5.0E+12	1.0E+13	7.8E+12
	3x F/T	ND	ND	ND	ND	ND	8.3E+12
	5x F/T	9.5E+12	ND	ND	ND	ND	8.3E+12
	10x F/T	9.9E+12	9.3E+09	700	6.2E+12	8.6E+12	7.8E+12
F1	Initial	8.3E+12	9.3E+09	1100	5.5E+12	9.5E+12	6.4E+12
	3x F/T	ND	ND	ND	ND	ND	6.7E+12
	5x F/T	1.1E+13	ND	ND	NO	ND	7.0E+12
	10x F/T	1.1E+13	2.9E+10	300	5.7E+12	9.5E+12	5.9E+12
F2	Initial	9.4E+12	1.1E+10	1100	5.2E+12	9.9E+12	7.3E+12
	3x F/T	ND	ND	ND	ND	ND	6.8E+12
	5x F/T	1.2E+13	ND	ND	ND	ND	1.0E+13
	10x F/T	1.2E+13	2.4E+10	400	5.6E+12	1.0E+13	6.4E+12
F3	Initial	7.7E+12	1.1E+09	8400	6.7E+12	9.7E+12	8.2E+12
	3x F/T	ND	ND	ND	ND	ND	6.1E+12
	5x F/T	1.8E+13	ND	ND	ND	ND	6.0E+12
	10x F/T	1.3E+13	3.6E+09	2600	5.7E+12	1.0E+13	5.7E+12

FIG. 4D




FIG. 5A




FIG. 5B

Formulation	Condition	Appearance	pH	Monomer Peak Mean Size (nm)	Polydisp. Index	% Monomer, DLS	% HMWS, DLS	% HMWS, SE-LC	% Full Capsids
F0	Initial	Pass	7.3	29.7	0.07	100.0%	0.0%	0.2%	89.2%
	1 week	Pass	7.3	ND	ND	ND	ND	ND	ND
	1 month	Pass	7.3	30.2	0.08	100.0%	0.0%	ND	ND
	2 months	Pass	7.2	30.1	0.07	100.0%	0.0%	0.7%	NAA
	3 months	Pass	7.4	32.9	0.13	100.0%	0.0%	0.0%	NAA
	6 months	Pass	7.2	30.6	0.08	100.0%	0.0%	ND	ND
F1	Initial	Pass	7.3	33.0	0.12	100.0%	0.0%	0.5%	89.1%
	1 week	Pass	7.1	ND	ND	ND	ND	ND	ND
	1 month	Pass	7.2	34.3	0.14	100.0%	0.0%	ND	ND
	2 months	Pass	7.2	34.4	0.14	100.0%	0.0%	1.9%	NAA
	3 months	Pass	7.2	36.4	0.17	99.6%	0.4%	5.5%	NAA
	6 months	Pass	7.2	ND	ND	ND	ND	ND	ND
F2	Initial	Pass	7.2	33.6	0.14	100.0%	0.0%	0.5%	89.3%
	1 week	Pass	7.2	ND	ND	ND	ND	ND	ND
	1 month	Pass	7.2	34.2	0.13	100.0%	0.0%	ND	ND
	2 months	Pass	7.2	34.6	0.14	100.0%	0.0%	0.8%	NAA
	3 months	Pass	7.2	35.4	0.17	99.5%	0.5%	1.3%	NAA
	6 months	Pass	7.3	34.0	0.13	100.0%	0.0%	ND	ND
F3	Initial	Pass	7.2	33.4	0.12	100.0%	0.0%	0.4%	89.1%
	1 week	Pass	7.2	ND	ND	ND	ND	ND	ND
	1 month	Pass	7.3	33.6	0.13	100.0%	0.0%	ND	ND
	2 months	Pass	7.2	33.5	0.12	100.0%	0.0%	0.6%	NAA
	3 months	Pass	7.4	34.7	0.17	99.2%	0.8%	1.7%	NAA
	6 months	Pass	7.2	33.9	0.13	100.0%	0.0%	ND	ND

FIG. 5C

Formulation	Condition	VG Titer (vg/mL)	Infectious Titer (TCID ₅₀ /mL)	VG/TCID ₅₀	Capsid Titer (capsids/mL)	Monomer Conc, SE-LC (p/mL)	Monomer Conc, MADLS (p/mL)
F0	Initial	7.9E+12	9.3E+09	1200	5.0E+12	1.0E+13	7.8E+12
	1 month	6.6E+12	ND	ND	ND	ND	8.6E+12
	2 months	4.9E+12	ND	ND	5.5E+12	9.0E+12	8.4E+12
	3 months	9.2E+12	2.4E+09	1100	5.3E+12	8.4E+12	8.6E+12
	6 months	ND	2.0E+09	ND	ND	ND	7.7E+12
F1	Initial	8.3E+12	9.3E+09	1100	5.5E+12	9.5E+12	6.4E+12
	1 month	4.5E+12	ND	ND	ND	ND	1.2E+13
	2 months	4.8E+12	ND	ND	5.3E+12	8.6E+12	1.5E+13
	3 months	1.1E+13	2.9E+09	600	6.3E+12	8.7E+12	4.8E+12
	6 months	ND	5.2E+08	ND	ND	ND	ND
F2	Initial	9.4E+12	1.1E+10	1100	5.2E+12	9.9E+12	7.3E+12
	1 month	3.9E+12	ND	ND	ND	ND	6.7E+12
	2 months	5.9E+12	ND	ND	5.5E+12	9.3E+12	7.6E+12
	3 months	8.2E+12	3.6E+09	1300	6.0E+12	1.1E+13	8.6E+12
	6 months	ND	2.4E+09	ND	ND	ND	6.5E+12
F3	Initial	7.7E+12	1.1 E+09	8400	6.7E+12	9.7E+12	8.2E+12
	1 month	4.8E+12	ND	ND	ND	ND	6.5E+12
	2 months	5.6E+12	ND	ND	5.3E+12	1.0E+13	7.1E +12
	3 months	9.5E+12	6.3E+09	500	6.8E+12	9.5E+12	7.3E+12
	6 months	ND	2.0E+09	ND	ND	ND	7.5E+12

FIG. 5D




FIG. 6A




FIG. 6B

FIG. 6C

Formulation	Condition	Appearance	pH	Monomer Peak Mean Size (nm)	Polydisp. Index	% Monomer, DLS	% HMWS, DLS	% HMWS, SE-LC	% Full Capsids
F0	Initial	Pass	7.3	29.7	0.07	100.0%	0.0%	0.2%	89.2%
	3 days	Pass	7.3	34.1	0.15	99.6 %	0.4%	ND	ND
	1 week	Pass	7.4	35.0	0.14	100.0%	0.0%	ND	ND
	2 weeks	Pass	7.3	35.9	0.20	98.1%	1.9%	ND	ND
	1 month	Pass	7.3	33.4	0.14	100.0%	0.0%	5.4%	NAA
	2 months	Pass	7.1	28.8	0.36	94.6%	5.4%	0.0%	NAA
	3 months	Fail	7.1	32.6	0.32	87.6%	12.4%	0.0%	NAA
	7 months	ND	ND	32.8	0.50	72.7%	27.3%	ND	ND
F1	Initial	Pass	7.3	33.0	0.12	100.0 %	0.0%	0.5%	89.1%
	3 days	Pass	7.2	ND	ND	ND	ND	ND	ND
	1 week	Pass	7.2	ND	ND	ND	ND	ND	ND
	2 weeks	Pass	7.2	ND	ND	ND	ND	ND	ND
	1 month	Pass	7.2	43.9	0.25	96.5%	3.6%	15.2%	NAA
	2 months	Pass	7.2	40.6	0.22	98.5%	1.5%	4.5%	NAA
	3 months	Fail	7.2	41.7	0.24	96.9%	3.1%	20.3%	NAA
	7 months	ND	ND	46.7	0.24	99.0%	1.0%	ND	ND
F2	Initial	Pass	7.2	33.6	0.14	100.0%	0.0%	0.5%	89.3%
	3 days	Pass	7.2	39.5	0.22	96.6%	3.4%	ND	ND
	1 week	Pass	7.2	45.1	0.25	94.5%	5.5%	ND	ND
	2 weeks	Pass	7.2	42.2	0.26	92.6%	7.4%	ND	ND
	1 month	Pass	7.2	41.3	0.23	97.6%	2.4%	6.5%	NAA
	2 months	Pass	7.2	42.0	0.23	98.5%	1.5%	11.6%	NAA
	3 months	Fail	7.3	40.3	0.24	95.0%	5.0%	17.1%	NAA
	7 months	ND	ND	55.4	0.28	97.5%	2.5%	ND	ND
F3	Initial	Pass	7.2	33.4	0.12	100.0%	0.0%	0.4%	89.1%
	3 days	Pass	7.2	36.3	0.23	96.6%	3.4%	ND	ND
	1 week	Pass	7.2	40.2	0.20	98.4%	1.6%	ND	ND
	2 weeks	Pass	7.2	42.9	0.22	98.2%	1.8%	ND	ND
	1 month	Pass	7.2	37.7	0.18	100.0%	0.0%	3.5%	NAA
	2 months	Pass	7.2	37.9	0.20	98.9%	1.1%	5.6%	NAA

Formulation	Condition	Appearance	pH	Monomer Peak Mean Size (nm)	Polydisp. Index	% Monomer, DLS	% HMWS, DLS	% HMWS, SE-LC	% Full Capsids
	3 months	Fail	7.3	39.2	0.25	92.8%	7.2%	9.5%	NAA
	7 months	ND	ND	46.9	0.36	82.4%	17.6%	ND	ND

FIG. 6C (cont'd)

FIG. 6D

Formulation	Condition	VG Titer (vg/mL)	Infectious Titer (TCID ₅₀ /mL)	VG/TCID ₅₀	Capsid Titer (capsids/mL)	Monomer Conc. SE-LC, (p/mL)	Monomer Conc. (MADLS (p/mL)
F0	Initial	7.9E+12	9.3E+09	1200	5.0E+12	1.0E+13	7.8 E+12
	3 days	ND	ND	ND	ND	ND	6.3E+12
	1 week	ND	ND	ND	ND	ND	5.4E+12
	2 weeks	ND	ND	ND	ND	ND	3.2E+12
	1 month	8.1E+12	1.7E+07	440000	4.5E+12	5.9E+12	4.0 E+12
	2 months	4.5E+12	ND	ND	3.2E+12	3.9E+12	5.1E+13
	3 months	ND	ND	ND	3.6E+12	4.4E+12	5.2E+12
F1	7 months	ND	ND	ND	ND	ND	1.2E+13
	Initial	8.3E+12	9.3E+09	1100	2.6E+12	9.5E+12	6.4E+12
	3 days	ND	ND	ND	ND	ND	ND
	1 week	ND	ND	ND	ND	ND	ND
	2 weeks	ND	ND	ND	ND	ND	ND
	1 month	3.8E+12	1.1E+07	280000	2.6E+12	4.2E+12	1.0E+13
	2 months	5.7E+12	ND	ND	2.5E+12	4.2E+12	8.6E+12
F2	3 months	ND	ND	ND	2.4E+12	6.6E+12	1.8E+13
	7 months	ND	ND	ND	ND	ND	1.0E+14
	Initial	9.4E+12	1.1E+10	1100	5.2E+12	9.9E+12	7.3E+12
	3 days	ND	ND	ND	ND	ND	9.2E+12
	1 week	ND	ND	ND	ND	ND	1.5E+13
	2 weeks	ND	ND	ND	ND	ND	1.5E+13
	1 month	7.8E+12	1.4E+07	410000	3.4E+12	4.9E+12	1.4E+13
F3	2 months	6.0E+12	ND	ND	2.9E+12	5.9E+12	2.2E+13
	3 months	ND	ND	ND	3.4E+12	9.0E+12	3.2E+13
	7 months	ND	ND	ND	ND	ND	1.0E+13
	Initial	7.7E+12	1.1E+09	8400	6.7E+12	9.7E+12	8.2E+12
	3 days	ND	ND	ND	ND	ND	6.8E+12
	1 week	ND	ND	ND	ND	ND	2.8E+13
	2 weeks	ND	ND	ND	ND	ND	1.5E+13
	1 month	8.7E+12	2.0E+08	42000	4.9E+12	6.2E+12	5.5E+12
	2 months	9.3E+12	ND	ND	2.7E+12	7.2E+12	8.0E+12

Formulation	Condition	VG Titer (vg/mL)	Infectious Titer (TCID ₅₀ /mL)	VG/TCID ₅₀	Capsid Titer (capsids/mL)	Monomer Conc. SE-LC, (p/mL)	Monomer Conc. (MADLS (p/mL)
	3 months	ND	ND	ND	2.8E+12	9.4E+12	1.3E+13
	7 months	ND	ND	ND	ND	ND	1.1E+13

FIG. 6D (cont'd)

Quality Attribute	Testing Method	Method Variability	Stability Acceptance Criteria	Failure Criteria
Appearance Liquid	Visual inspection	N/A	Clear colorless, free of particulates	Hazy, cloudy, or visual detection of particles
pH	Potentiometric	± 0.1	7.1 - 7.5	< 7.1 or > 7.5
Monomer Peak Mean Size	DLS	< 10%	< 10% change from initial (t_0)	> +25% change from initial
Polydispersity Index	DLS	< 10%	PDI < 0.20	PDI > 0.25
Monomer Fraction	DLS	< 2.5%	> 97.5% monomer	< 95% monomer
HMWS Fraction	DLS	< 2.5%	< 2.5% HMWS	> 5% HMWS
HMWS Fraction	SE-LC	< 10%	< 5% HMWS	> 10% HMWS
% Full Capsids	AEX-LC	< 10%	> 85% full capsids	< 80% full capsids
Vector Genome Titer	qPCR	36.0%	-36% to +36% change from initial	< -50% change from initial
Infectious Titer	TCID ₅₀	50.0%	-50% to +50% change from initial	< -60% change from initial
VG/TCID ₅₀ Ratio	Calculation	86.0%	-0.5 to +0.5 -log ₁₀ change from initial	< -1 -log ₁₀ change from initial
Capsid Concentration	ELISA	20.0%	-20% to +20% change from initial	< -40% change from initial
Monomer Concentration	SE-LC	30.0%	-30% to +30% change from initial	< -40% change from initial
Monomer Concentration	MADLS	< 30%	-25% to +25% change from initial	< -40% change from initial

FIG. 7

Form	Test Condition	Appearance, Liquid	pH	Monomer Peak Mean Size <i>% change from t₀</i>	Polydisp Index	Monomer Fraction, DLS	HMWS Fraction, DLS	HMWS Fraction, SE-LC	% Full Capsids, AEX-LC
F0	Initial	Pass/Fail	pH	% change from t ₀	PDI	% Monomer	% HMWS	% HMWS	% Full
	10x F/T	Pass	7.3	29.65	0.07	100.0%	0.0%	0.2%	89.2%
	3 mo 25°C / 60%RH	Fail	7.3	1.8%	0.07	99.6%	0.4%	0.0%	89.3%
	1 mo 40°C / 75%RH	Pass	7.4	11.1%	0.13	99.6%	0.4%	0.0%	NAA
	Assessment	Pass	7.3	12.6%	0.14	100.0%	0.0%	5.4%	NAA
F1	Initial	Fail	Pass	Neutral	Neutral	Pass	Pass	Neutral	Pass
	10x F/T	Pass	7.3	33.0	0.12	100.0%	0.0%	0.5%	89.1%
	3 mo 25°C / 60%RH	Pass	7.3	1.0%	0.12	100.0%	0.0%	0.6%	89.2%
	1 mo 40°C / 75%RH	Pass	7.2	10.3%	0.17	100.0%	0.0%	5.5%	NAA
	Assessment	Pass	7.2	32.9%	0.25	85.5%	14.5%	15.2%	NAA
F2	Initial	Pass	Pass	Fail	Fail	Fail	Fail	Fail	Pass
	10x F/T	Pass	7.2	33.6	0.14	100.0%	0.0%	0.5%	89.3%
	3 mo 25°C / 60%RH	Pass	7.2	0.0%	0.13	100.0%	0.0%	0.5%	89.3%
	1 mo 40°C / 75%RH	Pass	7.2	5.4%	0.17	98.4%	1.6%	1.3%	NAA
	Assessment	Pass	7.2	23.0%	0.23	99.7%	0.4%	6.5%	NAA
F3	Initial	Pass	Pass	Neutral	Neutral	Pass	Pass	Neutral	Pass
	10x F/T	Pass	7.2	33.4	0.12	100.0%	0.0%	0.4%	89.1%
	3 mo 25°C / 60%RH	Pass	7.3	1.0%	0.13	99.5%	0.5%	0.3%	89.0%
	1 mo 40°C / 75%RH	Pass	7.4	3.8%	0.17	100.0%	0.0%	1.7%	NAA
	Assessment	Pass	7.2	12.9%	0.18	100.0%	0.0%	3.5%	NAA
		Pass	Pass	Neutral	Neutral	Pass	Pass	Pass	Pass

FIG. 8

Form	Test	VG Titer, qPCR	Infect Titer, TCID50	VG / TCID50 Ratio	Capsid Conc, ELISA	Monomer Conc, SE-LC	Monomer Conc, MADLS
	Condition	% change from t0	% change from t0	Neg. log change from t0	% change from t0	% change from t0	% change from t0
F0	Initial	7.9E+12	9.3E+09	1200	5.0E+12	1.0E+13	7.8E+12
	10x F/T	25.1%	0.0%	0.23	24.5%	-13.3%	10.8%
	3 mo 25°C / 60%RH	16.9%	-73.9%	0.02	6.6%	-14.8%	32.0%
	1 mo 40°C / 75%RH	3.2%	-99.8%	-2.57	-9.6%	-40.7%	-43.3%*
	Assessment	Pass	No Score		Pass	Neutral	Neutral
F1	Initial	8.3E+12	9.3E+09	1100	5.5E+12	9.5E+12	6.4E+12
	10x F/T	36.5%	216.2%	0.54	4.2%	0.2%	-5.5%
	3 mo 25°C / 60%RH	31.4%	-68.3%	0.26	13.6%	-8.5%	-4.6%
	1 mo 40°C / 75%RH	-54.5%	-99.9%	-2.39	-52.1%	-55.5%	28.9%*
	Assessment	Fail	No Score		Neutral	Neutral	Neutral
F2	Initial	9.4E+12	1.1E+10	1100	5.2E+12	9.9E+12	7.3E+12
	10x F/T	27.3%	115.4%	0.46	7.7%	5.1%	-17.5%
	3 mo 25°C / 60%RH	-13.2%	-68.4%	-0.05	15.8%	6.8%	11.4%
	1 mo 40°C / 75%RH	-17.6%	-99.9%	-2.56	-34.4%	-50.7%	65.9%*
	Assessment	Pass	No Score		Pass	Neutral	Neutral
F3	Initial	7.7E+12	1.1E+09	8400	6.7E+12	9.7E+12	7.7E+12
	10x F/T	68.9%	216.2%	0.5	-15.2%	3.7%	-23.6%
	3 mo 25°C / 60%RH	22.8%	462.3%	1.27	2.6%	-1.5%	-15.0%
	1 mo 40°C / 75%RH	12.6%	-82.2%	-0.7	-27.0%	-36.1%	0.1%*
	Assessment	Neutral	No Score		Pass	Pass	Pass

* Result was not reproducible and not considered for final scoring

FIG. 9

Stability Results

Test	Temp (°C)	T0	1 month	3 months	6 months	9 months	12 months	18 months	24 months
Appearance, Liquid	-70	Clear	≤B9	≤B9	≤B9	≤B9	≤B9	≤B9	≤B9
	5		≤B9	≤B9	≤B9	≤B9	≤B9	≤B9	≤B9
	25		≤B9	≤B9	≤B9	≤B9	≤B9	≤B9	≤B9
	-70	Colorless	≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I
	5		≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I
	25		≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I	≤Ref I
	-70	No visible particulates	EFVP	EFVP	EFVP	EFVP	EFVP	EFVP	EFVP
	5		EFVP	EFVP	EFVP	EFVP	EFVP	EFVP	EFVP
	25		EFVP	EFVP	EFVP	EFVP	EFVP	EFVP	EFVP
Container Closure Integrity	-70							Pass	
pH	-70	7.4	7.5	7.4	7.4	7.4	7.4	7.3	7.3
	5		7.4	7.3	7.3	7.4	7.4		
	25		7.4	7.3					
Genome Titer (qPCR)	-70	1.06E+13	8.12E+12	9.17E+12	7.45E+12	1.14E+13	8.77E+12	9.05E+12	7.41E+12
	5		7.69E+12	8.93E+12	6.79E+12		9.00E+12		
	25		7.93E+12	9.26E+12					
Capsid Titer (ELISA)	-70	1.50E+13	1.05E+13	1.28E+13	1.07E+13	1.21E+13	1.40E+13	1.59E+13	1.07E+13
	5		1.14E+13	1.26E+13	1.14E+13		1.26E+13		
	25		1.06E+13	1.23E+13					
Infectious Titer (TCID ₅₀)	-70	4.31E+09	2.93E+09	2.20E+09	9.29E+08	3.56E+09	1.82E+09	2.00E+09	7.67E+09
	5		2.00E+09	9.27E+08	1.36E+09		5.22E+09		
	25		5.22E+09	1.12E+09					
Genome Titer / Infectious Titer Ratio	-70	2459	3618	4818	11410	2978	5824	5300	1382
Capsid Protein Purity (CE): % Capsid	-70	94.8	94.4	94.5	94.5	94.7	94.6	94.5	94.9
	5		94.2	94.8	94.5		94.8		
	25		94.3	93.5					
Capsid Protein Purity (CE): % Other	-70	NR	5.6	5.5	5.5	5.3	5.4	5.5	5.1
	5		5.8	5.2	5.5		5.2		
	25		5.7	6.5					

FIG. 10

Test	Temp (°C)	T0	1 month	3 months	6 months	9 months	12 months	18 months	24 months
Size Distribution (SEC): % Monomer	-70	99.4	98.9	98.9	99.2	99.0	99.3	99.2	99.0
	5		99.1	98.8	99.3		99.1		
	25		98.5	98.5					
Size Distribution (SEC): % Aggregates	-70	NMT 1.0	1.1	1.1	NMT 1.0	1.0	N/A	NMT 1.0	NMT 1.0
	5		NMT 1.0	1.2	NMT 1.0		NMT 1.0		
	25		1.5	1.5					
Particle Content (SEC): 260/280 Ratio	-70	1.32	1.32	1.32	1.32	1.32	1.32	1.33	1.32
	5		1.32	1.32	1.32		1.32		
	25		1.32	1.33					
In vitro FVIII Activity/ Functional Bioassay (% Relative Potency)	-70	73	76	82	64	60	69	60	55
	5		73	60	60		48		
	25		61	42					
P188 Concentration	-70	NT	443.8	454.4	456.8	432.3	464.2	504.8	399.0
	5		NT	467.2	455.7		466.6		
	25		454.7	461.9					

*1-5 particles (prctl); particle description – “1 long fiber particle”

NR: Not reported, NMT: Not more than, NT: not tested

FIG. 10 (cont'd)

International application No
PCT/US2021/064222

A. CLASSIFICATION OF SUBJECT MATTER INV. A61K9/08 A61K48/00 A61K47/26 A61P7/04 A61K47/02 A61K47/10 ADD. According to International Patent Classification (IPC) or to both national classification and IPC		
B. FIELDS SEARCHED Minimum documentation searched (classification system followed by classification symbols) A61K A61P Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched Electronic data base consulted during the international search (name of data base and, where practicable, search terms used) EPO-Internal, BIOSIS, EMBASE, WPI Data		
C. DOCUMENTS CONSIDERED TO BE RELEVANT		
Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	WO 2020/102369 A1 (REGENXBIO INC [US]) 22 May 2020 (2020-05-22) paragraphs [0223] - [0228], [0281] - [0282] -----	1-19, 22-33 20,21
A	WRIGHT J FRASER ET AL: "Recombinant adeno-associated virus: formulation challenges and strategies for a gene therapy vector", CURRENT OPINION IN DRUG DISCOVERY AND DEVELOPMENT, CURRENT DRUGS, LONDON, GB, vol. 6, no. 2, 1 March 2003 (2003-03-01), pages 174-178, XP009092903, ISSN: 1367-6733 abstract ----- -/--	1-33
☒ Further documents are listed in the continuation of Box C. ☒ See patent family annex.		
* Special categories of cited documents : "A" document defining the general state of the art which is not considered to be of particular relevance "E" earlier application or patent but published on or after the international filing date "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified) "O" document referring to an oral disclosure, use, exhibition or other means "P" document published prior to the international filing date but later than the priority date claimed "T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention "X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone "Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art "&" document member of the same patent family		
Date of the actual completion of the international search 23 March 2022		Date of mailing of the international search report 05/04/2022
Name and mailing address of the ISA/ European Patent Office, P.B. 5818 Patentlaan 2 NL - 2280 HV Rijswijk Tel. (+31-70) 340-2040, Fax: (+31-70) 340-3016		Authorized officer Konter, Jörg

INTERNATIONAL SEARCH REPORT

International application No

PCT/US2021/064222

C(Continuation). DOCUMENTS CONSIDERED TO BE RELEVANT		
Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	GNATENKO D V ET AL: "EXPRESSION OF THERAPEUTIC LEVELS OF FACTOR VIII IN HEMOPHILIA A MICE USING A NOVEL ADENO/ADENO-ASSOCIATED HYBRID VIRUS", THROMBOSIS AND HAEMOSTASIS, SCHATTAUER GMBH, DE, vol. 92, no. 2, 1 August 2004 (2004-08-01) , pages 317-327, XP009036414, ISSN: 0340-6245 abstract -----	1-33
X,P	WO 2021/071835 A1 (REGENXBIO INC [US]) 15 April 2021 (2021-04-15)	1-19, 22-33
A,P	paragraphs [0413] - [0414]; claim 1 -----	20,21

INTERNATIONAL SEARCH REPORT

International application No.

PCT/US2021/064222

Box No. I Nucleotide and/or amino acid sequence(s) (Continuation of item 1.c of the first sheet)

1. With regard to any nucleotide and/or amino acid sequence disclosed in the international application, the international search was carried out on the basis of a sequence listing:
 - a. ☒ forming part of the international application as filed:
 - ☒ in the form of an Annex C/ST.25 text file.
 - ☐ on paper or in the form of an image file.
 - b. ☐ furnished together with the international application under PCT Rule 13ter.1(a) for the purposes of international search only in the form of an Annex C/ST.25 text file.
 - c. ☐ furnished subsequent to the international filing date for the purposes of international search only:
 - ☐ in the form of an Annex C/ST.25 text file (Rule 13ter.1(a)).
 - ☐ on paper or in the form of an image file (Rule 13ter.1(b) and Administrative Instructions, Section 713).
2. ☐ In addition, in the case that more than one version or copy of a sequence listing has been filed or furnished, the required statements that the information in the subsequent or additional copies is identical to that forming part of the application as filed or does not go beyond the application as filed, as appropriate, were furnished.
3. Additional comments:

INTERNATIONAL SEARCH REPORT

Information on patent family members

International application No

PCT/US2021/064222

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 2020102369 A1	22-05-2020	AU 2019379141 A1 BR 112021009370 A2 CA 3117982 A1 EP 3880235 A1 JP 2022513034 A KR 20210092755 A SG 11202104295U A US 2021393802 A1 WO 2020102369 A1	03-06-2021 17-08-2021 22-05-2020 22-09-2021 07-02-2022 26-07-2021 29-06-2021 23-12-2021 22-05-2020
WO 2021071835 A1	15-04-2021	TW 202126319 A WO 2021071835 A1	16-07-2021 15-04-2021