

(19) United States

(12) Patent Application Publication (10) Pub. No.: US 2006/0199205 A1

Natsoulis et al.

Sep. 7, 2006 (43) Pub. Date:

(54) REAGENT SETS AND GENE SIGNATURES FOR RENAL TUBULE INJURY

(76) Inventors: Georges Natsoulis, Kensington, CA (US); Mark Fielden, Santa Clara, CA (US); Kurt Jarnagin, San Mateo, CA (US); Kyle Kolaja, San Mateo, CA (US)

Correspondence Address:

HOWREY LLP C/O IP DOCKETING DEPARTMENT 2941 FAIRVIEW PARK DRIVE, SUITE 200 **FALLS CHURCH, VA 22042-2924 (US)**

(21) Appl. No.: 11/357,887

(22) Filed: Feb. 17, 2006

Related U.S. Application Data

Continuation-in-part of application No. 11/184,272, filed on Jul. 18, 2005.

(60) Provisional application No. 60/589,409, filed on Jul. 19, 2004.

Publication Classification

(51) Int. Cl. (2006.01)C12Q 1/68 G01N 33/53 (2006.01)C12M 1/34 (2006.01)

(57)**ABSTRACT**

The invention discloses reagent sets and gene signatures for predicting onset of renal tubule injury in a subject. The invention also provides a necessary set of 186 genes useful for generating signatures of varying size and performance capable of predicting onset of renal tubule injury. The invention also provides methods, apparatuses and reagents useful for predicting future renal tubule injury based on expression levels of genes in the signatures. In one particular embodiment the invention provides a method for predict whether a compound will induce renal tubule injury using gene expression data from sub-acute treatments.

Accession	Weight (w) Average Ratio	Average Logio Ratio (1)	Average impact (w x r)	Percent Contribution	Unigene Identity
A105417	-0.89	-0.30	0.26	15.30	EST
BF 404557	. .38	-0.15	0.20	11.74	EST
U08257	0.88	0.17	0.15	8.73	Glutamate receptor, ionotropic, kainate 4
BF285022	1.46	0.10	0.15	8.61	EST
AF155910	0.55	0.23	0.12	7.22	heat shock 27kD protein family, member 7
A114646	0.83	0.17	0.11	6.13	EST, Similar to guanylate kinase
A105049	0.82	0.12	0.10	5.78	EST
4W916023	-0.64	0.11	20:0	4.21	EST, Similar to Kelch-like protein 8
A[227912	0.46	0.17	0.08	4.54	EST, Similar to Sorting nexin 3
3F403410	0.42	0.16	0.07	3.96	EST
700697	0.63	0.10	0.06	3.75	Cathepsin L
AW143082	-0.3	-0.17	90.0	2.94	EST
AL599126	0.36	0.12	0.04	2.52	EST, Similar to inner centromere protein-B
4102732	-0.31	-0.12	0.04	2.15	EST, Similar to microsomal signal peptidase 23 kDa subunit
41176933	0.46	89.0	0.04	2.04	EST
AF208288	-0.27	-0.13	0.04	2.06	G protein-coupled receptor 26
4F281635	0.43	0.05	0.02	1.29	zinc finger protein 22
J24174	0.09	0.19	0.02	0.99	cyclin-dependent kinase inhibitor 1A
AW142947	0 .22	0 .09	0.02	.	EST
BF396132	0.26	-0.05	0.01	0.78	echinoderm micratubule associated protein like 2
57049	-0.17	0 .09	0.01	0.85	methylenetetrahydrofolate reductase
UM_012610	90.0 0	-0.16	0.01	0.74	neive growth factor receptor
NW520754	0 .08	0 .1	0.01	0.53	potassium channel, subfamily K, member 3
AI231846	-0.13	9.0 .	0.01	0.46	EST
3E116947	0.05	0.13	0.01	0.36 0.36	EST
4W920B1B	0.03	0.23 6.23	0.0	0.35	EST
4W917933	0 .04	0 .13	0.01	0.3 8.0	EST
AW144517	0 .06	6 .0	0.005	0.28	adenylate cyclase 5
4B021980	-0.0 5	90.0	0.003	0.19	fatty acid desaturase 2
AF087454	O.29	-0.002	0.001	0.04	potassium voltage-gated channel, subfamily Q, member 3
3E097309	0.41	0.003	0.001	0.07	EST, Similar to BRPF1 protein
AW919837	0.05	0.01	-0.0005		EST
NM_13197	0.03	5. 25.	-0.01	•	aminolevulinic acid synthase 2
BF396955	0.72	0.07 20.00	0.06		EST, Similar to cell-cycle-dependent 350K nuclear protein
H-281149	ا ئ	-C.D2	-U.U.	1	EST, Similar to discs large nomotog /

Figure 1




Figure 1




Figure 2

REAGENT SETS AND GENE SIGNATURES FOR RENAL TUBULE INJURY

CROSS REFERENCE TO RELATED APPLICATIONS

[0001] This application is a continuation-in-part of U.S. patent application Ser. No. 11/184,272, filed on Jul. 18, 2005, which claims priority from U.S. Provisional Application No. 60/589,409, filed Jul. 19, 2004, each of which is hereby incorporated by reference herein in its entirety.

FIELD OF THE INVENTION

[0002] This invention relates to reagent sets and gene signatures useful for predicting the onset of renal tubule injury (RTI) in a subject. The invention also provides methods, apparatuses and kits useful for predicting occurrence of renal tubule injury based on expression levels of genes in the signatures. In one embodiment the invention provides a method for predicting whether a compound will induce renal tubule injury using gene expression data from sub-acute treatments.

BACKGROUND OF THE INVENTION

[0003] Renal tubule injury (also referred to herein as, "tubular nephrosis") is a common drug-induced toxicity that includes degenerative lesions of the renal tubules, such as acute tubular dilation, vacuolation and necrosis. Necrotic lesions of the tubules can arise as a consequence of septic, toxic or ischemic insult, and is a frequent cause of renal failure among hospitalized patients. Recognition is hampered by the lack of accurate markers and the shortcomings and over-reliance of serum markers of impaired glomerular filtration rate (i.e., serum creatinine and blood urea nitrogen) (see e.g., Schrier et al., "Acute renal failure: definitions, diagnosis, pathogenesis, and therapy," J Clin Invest, 114(1):5-14 (2004)). Drugs associated with the development of tubular nephrosis include aminoglycoside antibiotics, antifungals, antineoplastics, immunosuppresants and radiocontrast dyes, among others.

[0004] Similarly to the human clinical setting, long-term treatment of rats during preclinical drug development with relatively low doses of aminoglycoside antibiotics, heavy metal toxicants or antineoplastic drugs, for example, leads to the development of degenerative lesions of the renal tubules. However, histopathological or clinical indications of kidney injury are not readily apparent in the early course of treatment, thus necessitating expensive and lengthy studies.

[0005] The development of methods to predict the future onset of renal tubule injury (RTI) and gain a greater understanding of the underlying mechanism, would facilitate the development more reliable clinical diagnostics and safer therapeutic drugs. In addition, improved preclinical markers for RTI would dramatically reduce the time, cost, and amount of compound required in order to prioritize and select lead candidates for progression through drug development.

SUMMARY OF THE INVENTION

[0006] The present invention provides methods, reagent sets, gene sets, and associated apparatuses and kits, that allow one to determine the early onset of renal tubule injury (or nephrotoxicity) by measuring gene expression levels. In

one particular embodiment, the invention provides a RTI "necessary set" of 186 genes mined from a chemogenomic dataset. These genes are information-rich with respect to classifying biological samples for onset of RTI, even at sub-acute doses and time points of 5 days or earlier, where clinical and histopathological evidence of RTI are not manifested. Further, the invention discloses that the necessary set for RTI classification has the functional characteristic of reviving the performance of a fully depleted set of genes (for classifying RTI) by supplementation with random selections of as few as 10% of the genes from the set of 186. In addition, the invention discloses that selections from the necessary set made based on percentage impact of the selected genes may be used to generate high-performing linear classifiers for RTI that include as few as 4 genes. In one embodiment, the invention provides several different linear classifiers (or gene signatures) for RTI. For all of the disclosed embodiments based on the necessary set of 186 genes, the invention also provides reagent sets and kits comprising polynucleotides and/or polypeptides that represent a plurality of genes selected from the necessary set.

[0007] In one embodiment, the present invention provides a method for testing whether a compound will induce renal tubule injury in a test subject, the method comprising: administering a dose of a compound to at least one test subject; after a selected time period, obtaining a biological sample from the at least one test subject; measuring the expression levels in the biological sample of at least a plurality of genes selected from those listed in Table 4; determining whether the sample is in the positive class for renal tubule injury using a classifier comprising at least the plurality of genes for which the expression levels are measured. In one embodiment, the method is carried out wherein the test subject is a mammal selected from the group consisting of a human, cat, dog, monkey, mouse, pig, rabbit, and rat. In one preferred embodiment the test subject is a rat. In one embodiment, the biological sample comprises kidney tissue. In one embodiment, the method is carried out wherein the test compound is administered to the subject intravenously (IV), orally (PO, per os), or intraperitoneally (IP). In one embodiment, the method is carried out wherein the dose administered does not cause histological or clinical evidence of renal tubule injury at about 5 days, about 7 days, about 14 days, or even about 21 days. In one embodiment, the method is carried out wherein the expression levels are measured as log₁₀ ratios of compound-treated biological sample to a compound-untreated biological sample. In one embodiment, the method of the invention is carried out wherein the classifier is a linear classifier. In alternative embodiments, the classifier may be a non-linear classifier. In one embodiment, the method is carried out wherein the selected period of time is about 5 days or fewer, 7 days or fewer, 14 days or fewer, or even 21 days or fewer. In one embodiment of the method, the selected period of time is at least about 28 days.

[0008] In one embodiment, the method is carried out wherein the classifier comprises the genes and weights corresponding to any one of iterations 1 through 5 in Table 4. In one embodiment, the method of the invention is carried out wherein the classifier for renal tubule injury classifies each of the 64 compounds listed in Table 2 according to its label as nephrotoxic and non-nephrotoxic.

[0009] In one embodiment, the method is carried out wherein the linear classifier for renal tubule injury is capable of classifying a true label set with a log odds ratio at least 2 standard deviations greater than its performance classifying a random label set. In preferred embodiments of the method, the linear classifier for renal tubule injury is capable of performing with a training log odds ratio of greater than or equal to 4.35. In another embodiment, the plurality of genes includes at least 4 genes selected from those listed in Table 4, the four genes having at least having at least 2, 4, 8, 16, 32, or 64% of the total impact of all of the genes in Table 4.

[0010] The present invention also provides a gene sets, and reagent sets based on those gene sets, that are useful for testing whether renal tubule injury will occur in a test subject. In one embodiment, the invention provides a reagent set comprising a plurality of polynucleotides or polypeptides representing a plurality of genes selected from those listed in Table 4. In one embodiment, the reagent set comprises a plurality of genes includes at least 4 genes selected from those listed in Table 4, the 4 genes having at least 2% of the total impact of all of the genes in Table 4. In another embodiment, the reagent set comprises a plurality of genes includes at least 8 genes selected from those listed in Table 4, the 8 genes having at least 4% of the total impact of all of the genes in Table 4. Other embodiments include reagent sets based on subsets of genes randomly selected from Table 4, wherein the subset includes at least 4 genes having at least 1, 2, 4, 8, 16, 32, or 64% of the total impact. In preferred embodiments, the reagent sets of the invention include represent as few genes as possible from Table 4 while maximizing percentage of total impact. In preferred embodiments, the reagent sets of the invention include fewer than 1000, 500, 400, 300, 200, 100, 50, 20, 10, or even 8, polynucleotides or polypeptides representing the plurality of genes from Table 4. In one embodiment, the reagent sets consist essentially of polynucleotides or polypeptides representing the plurality of genes from Table 4. Further, the invention comprises kits comprising the reagent sets as components. In one embodiment, the reagent set is packaged in a single container consisting essentially of polynucleotides or polypeptides representing the plurality of genes from Table 4.

[0011] In one embodiment, the reagent sets of the invention comprise polynucleotides or polypeptides representing genes comprising a random selection of at least about 10% of the genes from Table 4, wherein the addition of said randomly selected genes to a fully depleted gene set for the renal tubule injury classification question increases the average logodds ratio of the linear classifiers generated by the depleted set to at least about 2.5. In another embodiment, a random selection of at least 20% of the genes from Table 4, wherein the addition of said randomly selected genes to a fully depleted gene set for the renal tubule injury classification question increases the average logodds ratio of the linear classifiers generated by the depleted set to at least about 3.3. In another embodiment, a random selection of at least 40% of the genes from Table 4, wherein the addition of said randomly selected genes to a fully depleted gene set for the renal tubule injury classification question increases the average logodds ratio of the linear classifiers generated by the depleted set to at least about 4.0. In other embodiments, reagent sets of the present invention comprise random selections of at least about 5%, 30%, 50%, 60%, 70%, 80%, 90%, or even 99% of the genes from Table 4, each which are capable of substantially increasing the average performance of a depleted set for generating classifiers RTI.

[0012] In one embodiment, the invention provides a reagent set for classifying renal tubule injury comprising a set of polynucleotides or polypeptides representing a plurality of genes selected from Table 4, wherein the addition of a random selection of at least 10% of said plurality of genes to the fully depleted set for the renal tubule injury classification question increases the average logodds ratio of the linear classifiers generated by the depleted set by at least 2-fold. In another embodiment, the reagent set includes at least 40% of said plurality of genes to the fully depleted set for the renal tubule injury classification question increases the average logodds ratio of the linear classifiers generated by the depleted set by at least 3-fold.

[0013] In another preferred embodiment the plurality of genes are selected from the variables of a linear classifier capable of classifying renal tubule injury with a training log odds ratio of greater than or equal to 4.35. In one preferred embodiment, the plurality of genes is the set of genes in any one of iterations 1 through 5 in Table 4. In another embodiment, the plurality of genes is the set of genes in any one of Tables 7, 8, 10, and 11. In one embodiment the reagents are polynucleotide probes capable of hybridizing to a plurality of genes selected from those listed in Table 4, and in a preferred embodiment, the polynucleotide probes are labeled.

[0014] In another embodiment, the reagents are primers for amplification of the plurality of genes. In one embodiment the reagents are polypeptides encoded by a plurality of genes selected from those listed in Table 4. Preferably the reagents are polypeptides that bind to a plurality proteins encoded by a plurality of genes selected from those listed in Table 4. In one preferred embodiment, the reagent set comprises secreted proteins encoded by genes listed in Table 4.

[0015] The present invention also provides an apparatus for predicting whether renal tubule injury will occur in a test subject comprising a reagent set as described above. In preferred embodiments, the apparatus comprises a device with reagents for detecting polynucleotides, wherein the reagents comprise or consist essentially of a reagent set for testing whether renal tubule injury will occur in a test subject as described above.

[0016] In one embodiment, the apparatus comprises at least a plurality of polynucleotides or polypeptides representing a plurality of genes selected from those listed in Table 4. In one embodiment the apparatus comprises a plurality of genes includes at least 4 genes selected from those listed in Table 4, the four genes having at least 2% of the total impact of the genes in Table 4. In another preferred embodiment the plurality of genes are variables in a linear classifier capable of classifying renal tubule injury with a training log odds ratio of greater than or equal to 4.35. In one embodiment, the apparatus comprises the plurality of genes listed in any one of iterations 1 through 5 in Table 4. In one preferred embodiment, the apparatus comprises polynucleotide probes capable of hybridizing to a plurality of genes selected from those listed in Table 4. In preferred embodiments, the apparatus comprises a plurality of polynucleotide probes bound to one or more solid surfaces. In one embodiment, the plurality of probes are bound to a single solid

surface in an array. Alternatively, the plurality of probes are bound to the solid surface on a plurality of beads. In another preferred embodiment, the apparatus comprises polypeptides encoded by a plurality of genes selected from those listed in Table 4. In one preferred embodiment, the polypeptides are secreted proteins encoded by genes listed in Table 4.

[0017] The present invention also provides a method for predicting renal tubule injury in an individual comprising: obtaining a biological sample from the individual after short-term treatment with compound; measuring the expression levels in the biological sample of at least a plurality of genes selected from Table 4; and determining whether the sample is in the positive class for renal tubule injury using a linear classifier comprising at least the plurality of genes for which the expression levels are measured; wherein a sample in the positive class indicates that the individual will have renal tubule injury following sub-chronic treatment with compound. In one preferred embodiment, the method for predicting renal tubule injury is carried out wherein the genes encode secreted proteins. In a preferred embodiment, the individual is a mammal, and preferably a rat. In another preferred embodiment, the biological sample is selected from blood, urine, hair or saliva. In another preferred embodiment of the method, the expression log₁₀ ratio is measured using an array of polynucleotides.

[0018] In another embodiment, the invention provides a method for monitoring treatment of an individual for renal tubule injury, or with a compound suspected of causing renal tubule injury, said method comprising: obtaining a biological sample from the individual after short-term treatment with compound; measuring the expression levels in the biological sample of at least a plurality of genes selected from Table 4; and determining whether the sample is in the positive class for renal tubule injury using a linear classifier comprising at least the plurality of genes for which the expression levels are measured; wherein a sample in the positive class indicates that the individual will have renal tubule injury. In a preferred embodiment, the individual is a mammal, and preferably a rat. In another preferred embodiment, the biological sample is selected from blood, urine, hair or saliva. In another preferred embodiment of the method, the expression \log_{10} ratio is measured using an array of polynucleotides.

BRIEF DESCRIPTION OF THE DRAWINGS

[0019] FIG. 1 depicts the 35 genes in the first iteration RTI signature derived according to the method of Example 3, their corresponding weights, and their average expression log, ratio in the 15 compound training positive class.

[0020] FIG. 2 depicts a plots of training and test logodds ratios for prediction of renal tubule injury for 20 subsets of genes randomly selected from the necessary set. A training or test LOR of 4.00 could be achieved by signatures of as few as 4 and 7 genes, respectively.

DETAILED DESCRIPTION OF THE INVENTION

I. Overview

[0021] The present invention provides methods for predicting whether compound treatments induce future renal

tubular injury following sub-chronic or long-term treatment using expression data from sub-acute or short-term treatments. The invention provides necessary and sufficient sets of genes and specific signatures comprising these genes that allow gene expression data to be used to identify the ability of a compound treatment to induce late onset renal tubule injury before the actual histological or clinical indication of the toxicity. Further, the invention provides reagent sets and diagnostic devices comprising the disclosed gene sets and signatures that may be used to deduce compound toxicity using short term studies, and avoiding lengthy and costly long term studies.

II. Definitions

[0022] "Multivariate dataset" as used herein, refers to any dataset comprising a plurality of different variables including but not limited to chemogenomic datasets comprising logratios from differential gene expression experiments, such as those carried out on polynucleotide microarrays, or multiple protein binding affinities measured using a protein chip. Other examples of multivariate data include assemblies of data from a plurality of standard toxicological or pharmacological assays (e.g., blood analytes measured using enzymatic assays, antibody based ELISA or other detection techniques).

[0023] "Variable" as used herein, refers to any value that may vary. For example, variables may include relative or absolute amounts of biological molecules, such as mRNA or proteins, or other biological metabolites. Variables may also include dosing amounts of test compounds.

[0024] "Classifier" as used herein, refers to a function of a set of variables that is capable of answering a classification question. A "classification question" may be of any type susceptible to yielding a yes or no answer (e.g., "Is the unknown a member of the class or does it belong with everything else outside the class?"). "Linear classifiers" refers to classifiers comprising a first order function of a set of variables, for example, a summation of a weighted set of gene expression logratios. A valid classifier is defined as a classifier capable of achieving a performance for its classification task at or above a selected threshold value. For example, a log odds ratio≥4.00 represents a preferred threshold of the present invention. Higher or lower threshold values may be selected depending of the specific classification task.

[0025] "Signature" as used herein, refers to a combination of variables, weighting factors, and other constants that provides a unique value or function capable of answering a classification question. A signature may include as few as one variable. Signatures include but are not limited to linear classifiers comprising sums of the product of gene expression logratios by weighting factors and a bias term.

[0026] "Weighting factor" (or "weight") as used herein, refers to a value used by an algorithm in combination with a variable in order to adjust the contribution of the variable.

[0027] "Impact factor" or "Impact" as used herein in the context of classifiers or signatures refers to the product of the weighting factor by the average value of the variable of interest. For example, where gene expression logratios are the variables, the product of the gene's weighting factor and the gene's measured expression \log_{10} ratio yields the gene's

impact. The sum of the impacts of all of the variables (e.g., genes) in a set yields the "total impact" for that set.

[0028] "Scalar product" (or "Signature score") as used herein refers to the sum of impacts for all genes in a signature less the bias for that signature. A positive scalar product for a sample indicates that it is positive for (i.e., a member of) the classification that is determined by the classifier or signature.

[0029] "Sufficient set" as used herein is a set of variables (e.g., genes, weights, bias factors) whose cross-validated performance for answering a specific classification question is greater than an arbitrary threshold (e.g., a log odds ratio ≥ 4.0).

[0030] "Necessary set" as used herein is a set of variables whose removal from the full set of all variables results in a depleted set whose performance for answering a specific classification question does not rise above an arbitrarily defined minimum level (e.g., log odds ratio ≥4.00).

[0031] "Log odds ratio" or "LOR" is used herein to summarize the performance of classifiers or signatures. LOR is defined generally as the natural log of the ratio of the odds of predicting a subject to be positive when it is positive, versus the odds of predicting a subject to be positive when it is negative. LOR is estimated herein using a set of training or test cross-validation partitions according to the following equation,

$$LOR = \ln \left(\sum_{i=1}^{c} TP_i + 0.5 \right) * \left(\sum_{i=1}^{c} TN_i + 0.5 \right)$$

$$\left(\sum_{i=1}^{c} FP_i + 0.5 \right) * \left(\sum_{i=1}^{c} FN_i + 0.5 \right)$$

where c (typically c=40 as described herein) equals the number of partitions, and TP_i , TN_i , FP_i , and FN_i represent the number of true positive, true negative, false positive, and false negative occurrences in the test cases of the i^{th} partition, respectively.

[0032] "Array" as used herein, refers to a set of different biological molecules (e.g., polynucleotides, peptides, carbohydrates, etc.). An array may be immobilized in or on one or more solid substrates (e.g., glass slides, beads, or gels) or may be a collection of different molecules in solution (e.g., a set of PCR primers). An array may include a plurality of biological polymers of a single class (e.g., polynucleotides) or a mixture of different classes of biopolymers (e.g., an array including both proteins and nucleic acids immobilized on a single substrate).

[0033] "Array data" as used herein refers to any set of constants and/or variables that may be observed, measured or otherwise derived from an experiment using an array, including but not limited to: fluorescence (or other signaling moiety) intensity ratios, binding affinities, hybridization stringency, temperature, buffer concentrations.

[0034] "Proteomic data" as used herein refers to any set of constants and/or variables that may be observed, measured or otherwise derived from an experiment involving a plurality of mRNA translation products (e.g., proteins, peptides, etc) and/or small molecular weight metabolites or exhaled gases associated with these translation products.

III. General Methods of the Invention

[0035] The present invention provides a method to derive multiple non-overlapping gene signatures for renal tubule injury. These non-overlapping signatures use different genes and thus each may be used independently in a predictive assay to confirm that an individual will suffer renal tubule injury. Furthermore, this method for identifying non-overlapping gene signatures also provides the list of all genes "necessary" to create a signature that performs above a certain minimal threshold level for a specific predicting renal tubule injury. This necessary set of genes also may be used to derive additional signatures with varying numbers of genes and levels of performance for particular applications (e.g., diagnostic assays and devices).

[0036] Classifiers comprising genes as variables and accompanying weighting factors may be used to classify large datasets compiled from DNA microarray experiments. Of particular preference are sparse linear classifiers. Sparse as used here means that the vast majority of the genes measured in the expression experiment have zero weight in the final linear classifier. Sparsity ensures that the sufficient and necessary gene lists produced by the methodology described herein are as short as possible. These short weighted gene lists (i.e., a gene signature) are capable of assigning an unknown compound treatment to one of two classes

[0037] The sparsity and linearity of the classifiers are important features. The linearity of the classifier facilitates the interpretation of the signature—the contribution of each gene to the classifier corresponds to the product of its weight and the value (i.e., log₁₀ ratio) from the microarray experiment. The property of sparsity ensures that the classifier uses only a few genes, which also helps in the interpretation. More importantly, the sparsity of the classifier may be reduced to a practical diagnostic apparatus or device comprising a relatively small set of reagents representing genes.

[0038] A. Gene Expression Related Datasets

[0039] a. Various Useful Data Types

[0040] The present invention may be used with a wide range of gene expression related data types to generate necessary and sufficient sets of genes useful for renal tubule injury signatures. In a preferred embodiment, the present invention utilizes data generated by high-throughput biological assays such as DNA microarray experiments, or proteomic assays. The datasets are not limited to gene expression related data but also may include any sort of molecular characterization information including, e.g., spectroscopic data (e.g., UV-Vis, NMR, IR, mass spectrometry, etc.), structural data (e.g., three-dimensional coordinates) and functional data (e.g., activity assays, binding assays). The gene sets and signatures produced by using the present invention may be applied in a multitude of analytical contexts, including the development and manufacture of detection devices (i.e., diagnostics).

[0041] b. Construction of a Gene Expression Dataset

[0042] The present invention may be used to identify necessary and sufficient sets of responsive genes within a gene expression dataset that are useful for predicting renal tubule injury. In a preferred embodiment, a chemogenomic dataset is used. For example, the data may correspond to

treatments of organisms (e.g., cells, worms, frogs, mice, rats, primates, or humans etc.) with chemical compounds at varying dosages and times followed by gene expression profiling of the organism's transcriptome (e.g., measuring mRNA levels) or proteome (e.g., measuring protein levels). In the case of multicellular organisms (e.g., mammals) the expression profiling may be carried out on various tissues of interest (e.g., liver, kidney, marrow, spleen, heart, brain, intestine). Typically, valid sufficient classifiers or signatures may be generated that answer questions relevant to classifying treatments in a single tissue type. The present specification describes examples of necessary and sufficient gene signatures useful for classifying chemogenomic data in liver tissue. The methods of the present invention may also be used however, to generate signatures in any tissue type. In some embodiments, classifiers or signatures may be useful in more than one tissue type. Indeed, a large chemogenomic dataset, like that exemplified in the present invention may reveal gene signatures in one tissue type (e.g., liver) that also classify pathologies in other tissues (e.g., intestine).

[0043] In addition to the expression profile data, the present invention may be useful with chemogenomic datasets including additional data types such as data from classic biochemistry assays carried out on the organisms and/or tissues of interest. Other data included in a large multivariate dataset may include histopathology, pharmacology assays, and structural data for the chemical compounds of interest.

[0044] One example of a chemogenomic multivariate dataset particularly useful with the present invention is a dataset based on DNA array expression profiling data as described in U.S. patent publication 2002/0174096 A1, published Nov. 21, 2002 (titled "Interactive Correlation of Compound Information and Genomic Information"), which is hereby incorporated by reference for all purposes. Microarrays are well known in the art and consist of a substrate to which probes that correspond in sequence to genes or gene products (e.g., cDNAs, mRNAs, cRNAs, polypeptides, and fragments thereof), can be specifically hybridized or bound at a known position. The microarray is an array (i.e., a matrix) in which each position represents a discrete binding site for a gene or gene product (e.g., a DNA or protein), and in which binding sites are present for many or all of the genes in an organism's genome.

[0045] As disclosed above, a treatment may include but is not limited to the exposure of a biological sample or organism (e.g., a rat) to a drug candidate (or other chemical compound), the introduction of an exogenous gene into a biological sample, the deletion of a gene from the biological sample, or changes in the culture conditions of the biological sample. Responsive to a treatment, a gene corresponding to a microarray site may, to varying degrees, be (a) upregulated, in which more mRNA corresponding to that gene may be present, (b) down-regulated, in which less mRNA corresponding to that gene may be present, or (c) unchanged. The amount of up-regulation or down-regulation for a particular matrix location is made capable of machine measurement using known methods (e.g., fluorescence intensity measurement). For example, a two-color fluorescence detection scheme is disclosed in U.S. Pat. Nos. 5,474,796 and 5,807,522, both of which are hereby incorporated by reference herein. Single color schemes are also well known in the art, wherein the amount of up- or down-regulation is determined in silico by calculating the ratio of the intensities from the test array divided by those from a control.

[0046] After treatment and appropriate processing of the microarray, the photon emissions are scanned into numerical form, and an image of the entire microarray is stored in the form of an image representation such as a color JPEG or TIFF format. The presence and degree of up-regulation or down-regulation of the gene at each microarray site represents, for the perturbation imposed on that site, the relevant output data for that experimental run or scan.

[0047] The methods for reducing datasets disclosed herein are broadly applicable to other gene and protein expression data. For example, in addition to microarray data, biological response data including gene expression level data generated from serial analysis of gene expression (SAGE, supra) (Velculescu et al., 1995, Science, 270:484) and related technologies are within the scope of the multivariate data suitable for analysis according to the method of the invention. Other methods of generating biological response signals suitable for the preferred embodiments include, but are not limited to: traditional Northern and Southern blot analysis; antibody studies; chemiluminescence studies based on reporter genes such as luciferase or green fluorescent protein; Lynx; READS (GeneLogic); and methods similar to those disclosed in U.S. Pat. No. 5,569,588 to Ashby et. al., "Methods for drug screening," the contents of which are hereby incorporated by reference into the present disclosure.

[0048] In another preferred embodiment, the large multivariate dataset may include genotyping (e.g., single-nucleotide polymorphism) data. The present invention may be used to generate necessary and sufficient sets of variables capable of classifying genotype information. These signatures would include specific high-impact SNPs that could be used in a genetic diagnostic or pharmacogenomic assay.

[0049] The method of generating classifiers from a multivariate dataset according to the present invention may be aided by the use of relational database systems (e.g., in a computing system) for storing and retrieving large amounts of data. The advent of high-speed wide area networks and the internet, together with the client/server based model of relational database management systems, is particularly well-suited for meaningfully analyzing large amounts of multivariate data given the appropriate hardware and software computing tools. Computerized analysis tools are particularly useful in experimental environments involving biological response signals (e.g., absolute or relative gene expression levels). Generally, multivariate data may be obtained and/or gathered using typical biological response signals. Responses to biological or environmental stimuli may be measured and analyzed in a large-scale fashion through computer-based scanning of the machine-readable signals, e.g., photons or electrical signals, into numerical matrices, and through the storage of the numerical data into relational databases. For example a large chemogenomic dataset may be constructed as described in U.S. patent publication 2005/0060102, published Mar. 17, 2005, which is hereby incorporated by reference for all purposes.

[0050] B. Generating Valid Gene Signatures from a Chemogenomic Dataset

[0051] a. Mining a Large Chemogenomic Dataset

[0052] Generally classifiers or signatures are generated (i.e., mined) from a large multivariate dataset by first labeling the full dataset according to known classifications and then applying an algorithm to the full dataset that produces a linear classifier for each particular classification question. Each signature so generated is then cross-validated using a standard split sample procedure.

[0053] The initial questions used to classify (i.e., the classification questions) a large multivariate dataset may be of any type susceptible to yielding a yes or no answer. The general form of such questions is: "Is the unknown a member of the class or does it belong with everything else outside the class?" For example, in the area of chemogenomic datasets, classification questions may include "modeof-action" questions such as "All treatments with drugs belonging to a particular structural class versus the rest of the treatments" or pathology questions such as "All treatments resulting in a measurable pathology versus all other treatments." In the specific case of chemogenomic datasets based on gene expression, it is preferred that the classification questions are further categorized based on the tissue source of the gene expression data. Similarly, it may be helpful to subdivide other types of large data sets so that specific classification questions are limited to particular subsets of data (e.g., data obtained at a certain time or dose of test compound). Typically, the significance of subdividing data within large datasets become apparent upon initial attempts to classify the complete dataset. A principal component analysis of the complete data set may be used to identify the subdivisions in a large dataset (see e.g., US 2003/0180808 A1, published Sep. 25, 2003, which is hereby incorporated by reference herein.) Methods of using classifiers to identify information rich genes in large chemogenomic datasets is also described in U.S. Ser. No. 11/114,998, filed Apr. 25, 2005, which is hereby incorporated by reference herein for all purposes.

[0054] Labels are assigned to each individual (e.g., each compound treatment) in the dataset according to a rigorous rule-based system. The +1 label indicates that a treatment falls in the class of interest, while a -1 label indicates that the variable is outside the class. Thus, with respect to the 64 compound treatments shown in Table 2 (see Example 2 below) used in generating an RTI signature, the "nephrotoxic" treatments were labeled +1, whereas the "non-nephrotoxic" were labeled -1. Information used in assigning labels to the various individuals to classify may include annotations from the literature related to the dataset (e.g., known information regarding the compounds used in the treatment), or experimental measurements on the exact same animals (e.g., results of clinical chemistry or histopathology assays performed on the same animal). A more detailed description of the general method for using classification questions to mine a chemogenomic dataset for signatures is described in U.S. Ser. No. 11/149,612, filed Jun. 10, 2005, and PCT/US2005/020695, filed Jun. 10, 2005, each of which is hereby incorporated in its entirety by reference herein.

[0055] b. Algorithms for Generating Valid Gene Signatures

[0056] Dataset classification may be carried out manually, that is by evaluating the dataset by eye and classifying the

data accordingly. However, because the dataset may involve tens of thousands (or more) individual variables, more typically, querying the full dataset with a classification question is carried out in a computer employing any of the well-known data classification algorithms.

[0057] In preferred embodiments, algorithms are used to query the full dataset that generate linear classifiers. In particularly preferred embodiments the algorithm is selected from the group consisting of: SPLP, SPLR and SPMPM. These algorithms are based respectively on Support Vector Machines (SVM), Logistic Regression (LR) and Minimax Probability Machine (MPM). They have been described in detail elsewhere (See e.g., El Ghaoui et al., op. cit; Brown, M. P., W. N. Grundy, D. Lin, N. Cristianini, C. W. Sugnet, T. S. Furey, M. Ares, Jr., and D. Haussler, "Knowledge-based analysis of microarray gene expression data by using support vector machines," *Proc Natl Acad Sci USA* 97: 262-267 (2000)).

[0058] Generally, the sparse classification methods SPLP, SPLR, SPMPM are linear classification algorithms in that they determine the optimal hyperplane separating a positive and a negative class. This hyperplane, H can be characterized by a vectorial parameter, w (the weight vector) and a scalar parameter, b (the bias): $H=\{x|w^Tx+b=0\}$.

[0059] For all proposed algorithms, determining the optimal hyperplane reduces to optimizing the error on the provided training data points, computed according to some loss function (e.g., the "Hinge loss," i.e., the loss function used in 1-norm SVMs; the "LR loss;" or the "MPM loss" augmented with a 1-norm regularization on the signature, w. Regularization helps to provide a sparse, short signature. Moreover, this 1-norm penalty on the signature will be weighted by the average standard error per gene. That is, genes that have been measured with more uncertainty will be less likely to get a high weight in the signature. Consequently, the proposed algorithms lead to sparse signatures, and take into account the average standard error information.

[0060] Mathematically, the algorithms can be described by the cost functions (shown below for SPLP, SPLR and SPMPM) that they actually minimize to determine the parameters w and b.

SPLP

$$\begin{split} \min_{w,b} \sum_{i} e_i + \rho \sum_{i} \sigma_i |w_i| s.t. y_i (w^T x_i + b) &\geq 1 - e_i \\ e_i &\geq 0, \, i = 1, \, \dots \, , \, N \end{split}$$

[0061] The first term minimizes the training set error, while the second term is the 1-norm penalty on the signature w, weighted by the average standard error information per gene given by sigma. The training set error is computed according to the so-called Hinge loss, as defined in the constraints. This loss function penalizes every data point that is closer than "1" to the separating hyperplane H, or is on the wrong side of H. Notice how the hyperparameter rho allows trade-off between training set error and sparsity of the signature w.

SPLR

$$\min_{w,b} \sum_{i} \log \left(1 + \exp\left(-y_i(w^Tx_i + b)\right)\right) + \rho \sum_{i} \sigma_i |w_i|$$

[0062] The first term expresses the negative log likelihood of the data (a smaller value indicating a better fit of the data), as usual in logistic regression, and the second term will give rise to a short signature, with rho determining the trade-off between both.

SPMPM

$$\min_{w} \sqrt{w^T \hat{\Gamma}_+ w} + \sqrt{w^T \hat{\Gamma}_- w} + \rho \sum_i \sigma_i |w_i| s.t. \ w^T (\hat{x}_+ - \hat{x}_-) = 1$$

[0063] Here, the first two terms, together with the constraint are related to the misclassification error, while the third term will induce sparsity, as before. The symbols with a hat are empirical estimates of the covariances and means of the positive and the negative class. Given those estimates, the misclassification error is controlled by determining w and b such that even for the worst-case distributions for the positive and negative class (which we do not exactly know here) with those means and covariances, the classifier will still perform well. More details on how this exactly relates to the previous cost function can be found in e.g., El Ghaoui, L., G. R. G. Lanckriet, and G. Natsoulis, 2003, "Robust classifiers with interval data" Report # UCB/CSD-03-1279. Computer Science Division (EECS), University of California, Berkeley, Calif.

[0064] As mentioned above, classification algorithms capable of producing linear classifiers are preferred for use with the present invention. In the context of chemogenomic datasets, linear classifiers may be used to generate one or more valid signatures capable of answering a classification question comprising a series of genes and associated weighting factors. Linear classification algorithms are particularly useful with DNA array or proteomic datasets because they provide simplified signatures useful for answering a wide variety of questions related to biological function and pharmacological/toxicological effects associated with genes or proteins. These signatures are particularly useful because they are easily incorporated into wide variety of DNA- or protein-based diagnostic assays (e.g., DNA microarrays).

[0065] However, some classes of non-linear classifiers, so called kernel methods, may also be used to develop short gene lists, weights and algorithms that may be used in diagnostic device development; while the preferred embodiment described here uses linear classification methods, it specifically contemplates that non-linear methods may also be suitable.

[0066] Classifications may also be carried using principle component analysis and/or discrimination metric algorithms well-known in the art (see e.g., US 2003/0180808 A1, published Sep. 25, 2003, which is hereby incorporated by reference herein).

[0067] Additional statistical techniques, or algorithms, are known in the art for generating classifiers. Some algorithms produce linear classifiers, which are convenient in many diagnostic applications because they may be represented as a weighted list of variables. In other cases non-linear classifier functions of the initial variables may be used. Other types of classifiers include decision trees and neural networks. Neural networks are universal approximators (Hornik, K., M. Stinchcombe, and H. White. 1989. "Multilayer feedforward networks are universal approximators, "Neural Networks 2: 359-366); they can approximate any measurable function arbitrarily well, and they can readily be used to model classification functions as well. They perform well on several biological problems, e.g., protein structure prediction, protein classification, and cancer classification using gene expression data (see, e.g., Bishop, C. M. 1996. Neural Networks for Pattern Recognition. Oxford University Press; Khan, J., J. S. Wei, M. Ringner, L. H. Saal, M. Ladanyi, F. Westermann, F. Berthold, M. Schwab, C. R. Antonescu, C. Peterson, and P. S. Meltzer. 2001. Classification and diagnostic prediction of cancers using gene expression profiling and artificial neural networks. Nat Med 7: 673-679; Wu, C. H., M. Berry, S. Shivakumar, and J. McLarty. 1995. Neural networks for full-scale protein sequence classification: sequence encoding with singular value decomposition. Machine Learning 21: 177-193).

[0068] c. Cross-Validation of Gene Signatures

[0069] Cross-validation of a gene signature's performance is an important step for determining whether the signature is sufficient. Cross-validation may be carried out by first randomly splitting the full dataset (e.g., a 60/40 split). A training signature is derived from the training set composed of 60% of the samples and used to classify both the training set and the remaining 40% of the data, referred to herein as the test set. In addition, a complete signature is derived using all the data. The performance of these signatures can be measured in terms of log odds ratio (LOR) or the error rate (ER) defined as:

 $LOR = \ln(((TP+0.5)*(TN+0.5))/((FP+0.5)*(FN+0.5)))$ and ER = (FP+FN)/N;

[0070] where TP, TN, FP, FN, and N are true positives, true negatives, false positives, false negatives, and total number of samples to classify, respectively, summed across all the cross validation trials. The performance measures are used to characterize the complete signature, the average of the training or the average of the test signatures.

[0071] The SVM algorithms described above are capable of generating a plurality of gene signatures with varying degrees of performance for the classification task. In order to identify that signatures that are to be considered "valid," a threshold performance is selected for the particular classification question. In one preferred embodiment, the classifier threshold performance is set as log odds ratio greater than or equal to 4.00 (i.e., LOR≥4.00). However, higher or lower thresholds may be used depending on the particular dataset and the desired properties of the signatures that are obtained. Of course many queries of a chemogenomic dataset with a classification question will not generate a valid gene signature.

[0072] Two or more valid gene signatures may be generated that are redundant or synonymous for a variety of reasons. Different classification questions (i.e., class definitions) may result in identical classes and therefore identical signatures. For instance, the following two class definitions define the exact same treatments in the database: (1) all treatments with molecules structurally related to statins; and (2) all treatments with molecules having an IC₅₀<1 μ M for inhibition of the enzyme HMG CoA reductase.

[0073] In addition, when a large dataset is queried with the same classification question using different algorithms (or even the same algorithm under slightly different conditions) different, valid signatures may be obtained. These different signatures may or may not comprise overlapping sets of variables; however, they each can accurately identify members of the class of interest.

[0074] For example, as illustrated in Table 1, two equally performing gene signatures (LOR=~7.0) for the fibrate class of compounds may be generated by querying a chemogenomic dataset with two different algorithms: SPLP and SPLR. Genes are designated by their accession number and a brief description. The weights associated with each gene are also indicated. Each signature was trained on the exact same 60% of the multivariate dataset and then cross validated on the exact same remaining 40% of the dataset. Both signatures were shown to exhibit the exact same level of performance as classifiers: two errors on the cross validation data set. The SPLP derived signature consists of 20 genes. The SPLR derived signature consists of eight genes. Only three of the genes from the SPLP signature are present in the eight gene SPLR signature.

[0075] Table 1: Two Gene Signatures for the Fibrate Class of Drugs

[0076] It is interesting to note that only three genes are common between these two signatures, (K03249, BF282712, and BF387347) and even those are associated with different weights. While many of the genes may be different, some commonalities may nevertheless be discerned. For example, one of the negatively weighted genes in the SPLP derived signature is NM_017136 encoding squalene epoxidase, a well-known cholesterol biosynthesis gene. Squalene epoxidase is not present in the SPLR derived signature but aceto-acteylCoA synthetase, another cholesterol biosynthesis gene is present and is also negatively weighted.

[0077] Additional variant signatures may be produced for the same classification task. For example, the average signature length (number of genes) produced by SPLP and SPLR, as well as the other algorithms, may be varied by use of the parameter p (see e.g., El Ghaoui, L., G. R. G. Lanckriet, and G. Natsoulis, 2003, "Robust classifiers with interval data" Report # UCB/CSD-03-1279. Computer Science Division (EECS), University of California, Berkeley, Calif.; and PCT publication WO 2005/017807 A2, published Feb. 24, 2005, each of which is hereby incorporated by reference herein). Varying p can produce signatures of different length with comparable test performance (Natsoulis et al., "Classification of a large microarray data set: Algorithm comparison and analysis of drug signatures, "Gen. Res. 15:724-736 (2005)). Those signatures are obviously different and often have no common genes between them (i.e., they do not overlap in terms of genes used).

 $\cite{[0078]}$ C. "Stripping" Signatures from a Dataset to Generate the "Necessary" Set

[0079] Each individual classifier or signature is capable of classifying a dataset into one of two categories or classes

	Accession	Weight Unigene name
RLPC	K03249	1.1572 enoyl-Co A, hydratase/3-hydroxyacyl Co A dehydrogenase
	AW916833	1.0876 hypothetical protein RMT-7
	BF387347	0.4769 ESTs
	BF282712	0.4634 ESTs
	AF034577	0.3684 pyruvate dehydrogenate kinase 4
	NM_019292	0.3107 carbonic anhydrase 3
	AI179988	0.2735 ectodermal-neural cortex (with BTB-like domain)
	AI715955	0.211 Stac protein (SRC homology 3 and cysteine-rich domain protein)
	BE110695	0.2026 activating transcription factor 1
	J03752	0.0953 microsomal glutathione S-transferase 1
	D86580	0.0731 nuclear receptor subfamily 0, group B, member 2
	BF550426	0.0391 KDEL (Lys-Asp-Glu-Leu) endoplasmic reticulum protein retention receptor 2
	AA818999	0.0296 muscleblind-like 2
	NM_019125	0.0167 probasin
	AF150082	-0.0141 translocase of inner mitochondrial membrane 8 (yeast) homolog A
	BE118425	-0.0781 Arsenical pump-driving ATPase
	NM_017136	-0.126 squalene epoxidase
	AI171367	-0.3222 HSPC154 protein
	NM_019369	-0.637 inter alpha-trypsin inhibitor, heavy chain 4
	AI137259	-0.7962 ESTs
SPLR	NM_017340	•
	BF282712	4.1052 ESTs
	NM_012489	, ,
	BF387347	1.767 ESTs
	K03249	1.7524 enoyl-Co A, hydratase/3-hydroxyacyl Co A dehydrogenase
	NM_016986	, , , ,
	AB026291	-0.7456 acetoacetyl-CoA synthetase
	AI454943	-1.6738 likely ortholog of mouse porcupine homolog

defined by the classification question. Typically, an individual signature with the highest test log odds ratio will be considered as the best classifier for a given task. However, often the second, third (or lower) ranking signatures, in terms of performance, may be useful for confirming the classification of compound treatment, especially where the unknown compound yields a borderline answer based on the best classifier. Furthermore, the additional signatures may identify alternative sources of informational rich data associated with the specific classification question. For example, a slightly lower ranking gene signature from a chemogenomic dataset may include those genes associated with a secondary metabolic pathway affected by the compound treatment. Consequently, for purposes of fully characterizing a class and answering difficult classification questions, it is useful to define the entire set of variables that may be used to produce the plurality of different classifiers capable of answering a given classification question. This set of variables is referred to herein as a "necessary set." Conversely, the remaining variables from the full dataset are those that collectively cannot be used to produce a valid classifier, and therefore are referred to herein as the "depleted set."

[0080] The general method for identifying a necessary set of variables useful for a classification question involved what is referred to herein as a classifier "stripping" algorithm. The stripping algorithm comprises the following steps: (1) querying the full dataset with a classification question so as to generate a first linear classifier capable of performing with a log odds ratio greater than or equal to 4.0 comprising a first set of variables; (2) removing the variables of the first linear classifier from the full dataset thereby generating a partially depleted dataset; (3) re-querying the partially depleted dataset with the same classification question so as to generate a second linear classifier and crossvalidating this second classifier to determine whether it performs with a log odds ratio greater than or equal to 4. If it does not, the process stops and the dataset is fully depleted for variables capable of generating a classifier with an average log odds ratio greater than or equal to 4.0. If the second classifier is validated as performing with a log odds ratio greater than or equal to 4.0, then its variables are stripped from the full dataset and the partially depleted set if re-queried with the classification question. These cycles of stripping and re-querying are repeated until the performance of any remaining set of variables drops below an arbitrarily set LOR. The threshold at which the iterative process is stopped may be arbitrarily adjusted by the user depending on the desired outcome. For example, a user may choose a threshold of LOR=0. This is the value expected by chance alone. Consequently, after repeated stripping until LOR=0 there is no classification information remaining in the depleted set. Of course, selecting a lower value for the threshold will result in a larger necessary set.

[0081] Although a preferred cut-off for stripping classifiers is LOR=4.0, this threshold is arbitrary. Other embodiments within the scope of the invention may utilize higher or lower stripping cutoffs e.g., depending on the size or type of dataset, or the classification question being asked. In addition other metrics could be used to assess the performance (e.g., specificity, sensitivity, and others). Also the stripping algorithm removes all variables from a signature if it meets the cutoff. Other procedures may be used within the scope of the invention wherein only the highest weighted or ranking variables are stripped. Such an approach based on

variable impact would likely result in a classifier "surviving" more cycles and defining a smaller necessary set.

[0082] Other procedures may be used within the scope of the invention wherein only the highest weighted or ranking variables are stripped. Such an approach based on variable impact would likely result in a classifier "surviving" more cycles and defining a smaller necessary set.

[0083] In another alternative approach, the genes from signatures may be stripped from the dataset until it is unable to generate a signature capable of classifying the "true label set" with an LOR that is statistically different from its classification of the "random label set." The "true label set" refers to a training set of compound treatment data that is correctly labeled (e.g., +1 class, -1 class) for the particular classification question. The "random label set" refers to the same set of compound treatment data where the class labels have been randomly assigned. Attempts to use a signature to classify a random label set will result in an average LOR of approximately zero and some standard deviation (SD). These values may be compared to the average LOR and SD for the classifying the true label set, where the SD is calculated based on LOR results across the 20 or 40 splits. The difference in classifying true and random label sets with valid signatures should be significantly greater than random. In such an alternative approach, the selected performance threshold for a signature is a p-value rather than a LOR cutoff.

[0084] The resulting fully-depleted set of variables that remains after a classifier is fully stripped from the full dataset cannot generate a classifier for the specific classification question (with the desired level of performance). Consequently, the set of all of the variables in the classifiers that were stripped from the full set are defined as "necessary" for generating a valid classifier.

[0085] The stripping method utilizes a classification algorithm at its core. The examples presented here use SPLP for this task. Other algorithms, provided that they are sparse with respect to genes could be employed. SPLR and SPMPM are two alternatives for this functionality (see e.g., El Ghaoui, L., G. R. G. Lanckriet, and G. Natsoulis, 2003, "Robust classifiers with interval data" *Report # UCB/CSD*-03-1279. Computer Science Division (EECS), University of California, Berkeley, Calif., and PCT publication WO 2005/017807 A2, published Feb. 24, 2005, which is hereby incorporated by reference herein).

[0086] In one embodiment, the stripping algorithm may be used on a chemogenomics dataset comprising DNA microarray data. The resulting necessary set of genes comprises a subset of highly informative genes for a particular classification question. Consequently, these genes may be incorporated in diagnostic devices (e.g., polynucleotide arrays) where that particular classification (e.g., renal tubule injury) is of interest. In other exemplary embodiments, the stripping method may be used with datasets from proteomic experiments

[0087] D. Mining the Renal Tubule Injury Necessary Set for Signatures

[0088] Besides identifying the "necessary" set of genes for a particular signature (i.e., classifier), another important use of the stripping algorithm is the identification of multiple, non-overlapping sufficient sets of genes useful for answering

a particular classification question. These non-overlapping sufficient sets are a direct product of the above-described general method of stripping valid classifiers. Where the application of the method results in a second validated classifier with the desired level of performance, that second classifier by definition does not include any genes in common with the first classifier. Typically, the earlier stripped non-overlapping gene signature yields higher performance with fewer genes. In other words, the earliest identified sufficient set usually comprises the highest impact, most information-rich genes with respect to the particular classification question. The valid classifiers that appear during later iterations of the stripping algorithm typically contain a larger number of genes. However, these later appearing classifiers may provide valuable information regarding normally unrecognized relationships between genes in the dataset. For example, in the case of non-overlapping gene signatures identified by stripping in a chemogenomics dataset, the later appearing signatures may include families of genes not previously recognized as involved in the particular metabolic pathway that is being affected by a particular compound treatment. Thus, functional analysis of a gene signature stripping procedure may identify new metabolic targets associated with a compound treatment.

[0089] The necessary set high impact genes generated by the stripping method itself represents a subset of genes that may be mined for further signatures. Hence, the complete set of genes in a necessary set for predicting renal tubule injury may used to randomly generate random subsets of genes of varying size that are capable of generating additional predictive signatures. One preferred method of selecting such subsets is based on percentage of total impact. Thus, subsets of genes are selected whose summed impact factors are a selected percentage of the total impact (i.e., the sum of the impacts of all genes in the necessary set). These percentage impact subsets may be used to generate new signatures for predicting renal tubule injury. For example, a random subset from the necessary set of 9 genes with 4% of the total impact may be used with one of the SVM algorithms to generate a new linear classifier of 8 genes, weighting factors and a bias term that may be used as a signature for renal tubule injury. Thus, the necessary set for a particular classification represents a greatly reduced dataset that can generate new signatures with varying properties such as shorter (or longer) gene lengths and higher (or lower) LOR performance values.

[0090] E. Functional Characterization of the Renal Tubule Injury Necessary Set

[0091] The stripping method described herein produces a necessary set of genes representing for answering the RTI classification question. The RTI necessary set of genes also may be characterized in functional terms based on the ability of the information rich genes in the set to supplement (i.e., "revive") the ability of a fully "depleted" set of genes to generate valid RTI signatures. Thus, the necessary set for the RTI classification question corresponds to that set of genes from which any random selection when added to a depleted set (i.e., depleted for RTI classification question) restores the ability of that set to produce RTI signatures with an average LOR (avg. LOR) above a threshold level. The general method for functionally characterizing a necessary set in terms of its ability to revive its depleted set is described in U.S. Ser. No. 11/149,612, filed Jun. 10, 2005, and PCT/

US2005/020695, filed Jun. 10, 2005, each of which is hereby incorporated in its entirety by reference herein.

[0092] Preferably, the threshold performance used is an avg. LOR greater than or equal to 4.00. Other values for performance, however, may be set. For example, avg. LOR may vary from about 1.0 to as high as 8.0. In preferred embodiments, the avg. LOR threshold may be 3.0 to as high as 7.0 including all integer and half-integer values in that range. The necessary set may then be defined in terms of percentage of randomly selected genes from the necessary set that restore the performance of a depleted set above a certain threshold. Typically, the avg. LOR of the depleted set is ~1.20, although as mentioned above, datasets may be depleted more or less depending on the threshold set, and depleted sets with avg. LOR as low as 0.0 may be used. Generally, the depleted set will exhibit an avg. LOR between about 0.5 and 1.5.

[0093] The third parameter establishing the functional characteristics of the RTI necessary set of genes for answering the RTI classification question is the percentage of randomly selected genes from that set that result in reviving the threshold performance of the depleted set. Typically, where the threshold avg. LOR is at least 4.00 and the depleted set performs with an avg. LOR of ~1.20, typically 16-36% of randomly selected genes from the necessary set are required to restore the average performance of the depleted set to the threshold value. In preferred embodiments, the random supplementation may be achieved using 16, 18, 20, 22, 24, 26, 28, 30, 32, 34 or 36% of the necessary set.

[0094] Alternatively, as described above, the necessary set may be characterized based on its ability to randomly generate signatures capable of classifying a true label set with an average performance above those signatures ability to classify a random label set. In preferred embodiments, signatures generated from a random selection of at least 10% of the genes in the necessary set may perform at least 1 standard deviation, and preferably at least 2 standard deviations, better for classifying the true versus the random label set. In other embodiments, the random selection may be of at least 15%, 20%, 25%, 30%, 40%, 50%, and even higher percentages of genes from the set.

[0095] F. Using Signatures and the Necessary Set to Generate Diagnostic Assays and Devices for Predicting Renal Tubule Injury

[0096] A diagnostic usually consists in performing one or more assays and in assigning a sample to one or more categories based on the results of the assay(s). Desirable attributes of a diagnostic assays include high sensitivity and specificity measured in terms of low false negative and false positive rates and overall accuracy. Because diagnostic assays are often used to assign large number of samples to given categories, the issues of cost per assay and throughput (number of assays per unit time or per worker hour) are of paramount importance.

[0097] Typically the development of a diagnostic assay involves the following steps: (1) define the end point to diagnose, e.g., cholestasis, a pathology of the liver (2) identify one or more markers whose alteration correlates with the end point, e.g., elevation of bilirubin in the blood-stream as an indication of cholestasis; and (3) develop a

specific, accurate, high-throughput and cost-effective assay for that marker. In order to increase throughput and decrease costs several diagnostics are often combined in a panel of assays, especially when the detection methodologies are compatible. For example several ELISA-based assays, each using different antibodies to ascertain different end points may be combined in a single panel and commercialized as a single kit. Even in this case, however, each of the ELISA-based assays had to be developed individually often requiring the generation of specific reagents.

[0098] The present invention provides signatures and methods for identifying additional signatures comprising as few as 4 genes that are useful for determining a therapeutic or toxicological end-point for renal tubule injury. These signatures (and the genes from which they are composed) may also be used in the design of improved diagnostic devices that answer the same questions as a large microarray but using a much smaller fraction of data. Generally, the reduction of information in a large chemogenomic dataset to a simple signature enables much simpler devices compatible with low cost high throughput multi-analyte measurement.

[0099] As described herein, a large chemogenomic dataset may be mined for a plurality of informative genes useful for answering classification questions. The size of the classifiers or signatures so generated may be varied according to experimental needs. In addition, multiple non-overlapping classifiers may be generated where independent experimental measures are required to confirm a classification. Generally, the sufficient classifiers result in a substantial reduction of data that needs to be measured to classify a sample. Consequently, the signatures and methods of the present invention provide the ability to produce cheaper, higher throughput, diagnostic measurement methods or strategies. In particular, the invention provides diagnostic reagent sets useful in diagnostic assays and the associated diagnostic devices and kits. As used herein, diagnostic assays includes assays that may be used for patient prognosis and therapeutic monitoring.

[0100] Diagnostic reagent sets may include reagents representing the subset of genes found in the necessary set of 186 consisting of less than 50%, 40%, 30%, 20%, 10%, or even less than 5% of the total genes. In one preferred embodiment, the diagnostic reagent set is a plurality of polynucleotides or polypeptides representing specific genes in a sufficient or necessary set of the invention. Such biopolymer reagent sets are immediately applicable in any of the diagnostic assay methods (and the associate kits) well known for polynucleotides and polypeptides (e.g., DNA arrays, RT-PCR, immunoassays or other receptor based assays for polypeptides or proteins). For example, by selecting only those genes found in a smaller yet "sufficient" gene signature, a faster, simpler and cheaper DNA array may be fabricated for that signature's specific classification task. Thus, a very simple diagnostic array may be designed that answers 3 or 4 specific classification questions and includes only 60-80 polynucleotides representing the approximately 20 genes in each of the signatures. Of course, depending on the level of accuracy required the LOR threshold for selecting a sufficient gene signature may be varied. A DNA array may be designed with many more genes per signature if the LOR threshold is set at e.g., 7.00 for a given classification question. The present invention includes diagnostic devices based on gene signatures exhibiting levels of performance varying from less than LOR=3.00 up to LOR=10.00 and greater.

[0101] The diagnostic reagent sets of the invention may be provided in kits, wherein the kits may or may not comprise additional reagents or components necessary for the particular diagnostic application in which the reagent set is to be employed. Thus, for a polynucleotide array applications, the diagnostic reagent sets may be provided in a kit which further comprises one or more of the additional requisite reagents for amplifying and/or labeling a microarray probe or target (e.g., polymerases, labeled nucleotides, and the like).

[0102] A variety of array formats (for either polynucleotides and/or polypeptides) are well-known in the art and may be used with the methods and subsets produced by the present invention. In one preferred embodiment, photolithographic or micromirror methods may be used to spatially direct light-induced chemical modifications of spacer units or functional groups resulting in attachment at specific localized regions on the surface of the substrate. Light-directed methods of controlling reactivity and immobilizing chemical compounds on solid substrates are well-known in the art and described in U.S. Pat. Nos. 4,562,157, 5,143,854, 5,556,961, 5,968,740, and 6,153,744, and PCT publication WO 99/42813, each of which is hereby incorporated by reference herein.

[0103] Alternatively, a plurality of molecules may be attached to a single substrate by precise deposition of chemical reagents. For example, methods for achieving high spatial resolution in depositing small volumes of a liquid reagent on a solid substrate are disclosed in U.S. Pat. Nos. 5,474,796 and 5,807,522, both of which are hereby incorporated by reference herein.

[0104] It should also be noted that in many cases a single diagnostic device may not satisfy all needs. However, even for an initial exploratory investigation (e.g., classifying drug-treated rats) DNA arrays with sufficient gene sets of varying size (number of genes), each adapted to a specific follow-up technology, can be created. In addition, in the case of drug-treated rats, different arrays may be defined for each tissue

[0105] Alternatively, a single substrate may be produced with several different small arrays of genes in different areas on the surface of the substrate. Each of these different arrays may represent a sufficient set of genes for the same classification question but with a different optimal gene signature for each different tissue. Thus, a single array could be used for particular diagnostic question regardless of the tissue source of the sample (or even if the sample was from a mixture of tissue sources, e.g., in a forensic sample).

[0106] In addition, it may be desirable to investigate classification questions of a different nature in the same tissue using several arrays featuring different non-overlapping gene signatures for a particular classification question.

[0107] As described above, the methodology described here is not limited to chemogenomic datasets and DNA microarray data. The invention may be applied to other types of datasets to produce necessary and sufficient sets of variables useful for classifiers. For example, proteomics assay techniques, where protein levels are measured or

protein interaction techniques such as yeast 2-hybrid or mass spectrometry also result in large, highly multivariate dataset, which could be classified in the same way described here. The result of all the classification tasks could be submitted to the same methods of signature generation and/or classifier stripping in order to define specific sets of proteins useful as signatures for specific classification questions.

[0108] In addition, the invention is useful for many traditional lower throughput diagnostic applications. Indeed the invention teaches methods for generating valid, high-performance classifiers consisting of 5% or less of the total variables in a dataset. This data reduction is critical to providing a useful analytical device. For example, a large chemogenomic dataset may be reduced to a signature comprising less than 5% of the genes in the full dataset. Further reductions of these genes may be made by identifying only those genes whose product is a secreted protein. These secreted proteins may be identified based on known annotation information regarding the genes in the subset. Because the secreted proteins are identified in the sufficient set useful as a signature for a particular classification question, they are most useful in protein based diagnostic assays related to that classification. For example, an antibody-based blood serum assay may be produced using the subset of the secreted proteins found in the sufficient signature set. Hence, the present invention may be used to generate improved protein-based diagnostic assays from DNA array informa-

[0109] The general method of the invention as described above is exemplified below. The following examples are offered as illustrations of specific embodiments and are not intended to limit the inventions disclosed throughout the whole of the specification.

EXAMPLES

Example 1

Construction of Chemogenomic Reference Database (DrugMatrixTM)

[0110] This example illustrates the construction of a large multivariate chemogenomic dataset based on DNA microarray analysis of rat tissues from over 580 different in vivo compound treatments. This dataset was used to generate RTI signatures comprising genes and weights which subsequently were used to generate a necessary set of highly responsive genes that may be incorporated into high throughput diagnostic devices as described in Examples 2-7.

[0111] The detailed description of the construction of this chemogenomic dataset is described in Examples 1 and 2 of Published U.S. Pat. Appl. No. 2005/0060102 A1, published Mar. 17, 2005, which is hereby incorporated by reference for all purposes. Briefly, in vivo short-term repeat dose rat studies were conducted on over 580 test compounds, including marketed and withdrawn drugs, environmental and industrial toxicants, and standard biochemical reagents. Rats (three per group) were dosed daily at either a low or high dose. The low dose was an efficacious dose estimated from the literature and the high dose was an empirically-determined maximum tolerated dose, defined as the dose that causes a 50% decrease in body weight gain relative to controls during the course of the 5 day range finding study.

Animals were necropsied on days 0.25, 1, 3, and 5 or 7. Up to 13 tissues (e.g., liver, kidney, heart, bone marrow, blood, spleen, brain, intestine, glandular and nonglandular stomach, lung, muscle, and gonads) were collected for histopathological evaluation and microarray expression profiling on the Amersham CodeLink™ RU1 platform. In addition, a clinical pathology panel consisting of 37 clinical chemistry and hematology parameters was generated from blood samples collected on days 3 and 5.

[0112] In order to assure that all of the dataset is of high quality a number of quality metrics and tests are employed. Failure on any test results in rejection of the array and exclusion from the data set. The first tests measure global array parameters: (1) average normalized signal to background, (2) median signal to threshold, (3) fraction of elements with below background signals, and (4) number of empty spots. The second battery of tests examines the array visually for unevenness and agreement of the signals to a tissue specific reference standard formed from a number of historical untreated animal control arrays (correlation coefficient>0.8). Arrays that pass all of these checks are further assessed using principle component analysis versus a dataset containing seven different tissue types; arrays not closely clustering with their appropriate tissue cloud are discarded.

[0113] Data collected from the scanner is processed by the Dewarping/DetrendingTM normalization technique, which uses a non-linear centralization normalization procedure (see, Zien, A., T. Aigner, R. Zimmer, and T. Lengauer. 2001. Centralization: A new method for the normalization of gene expression data. *Bioinformatics*) adapted specifically for the CodeLink microarray platform. The procedure utilizes detrending and dewarping algorithms to adjust for non-biological trends and non-linear patterns in signal response, leading to significant improvements in array data quality.

[0114] Log₁₀-ratios are computed for each gene as the difference of the averaged logs of the experimental signals from (usually) three drug-treated animals and the averaged logs of the control signals from (usually) 20 mock vehicletreated animals. To assign a significance level to each gene expression change, the standard error for the measured change between the experiments and controls is computed. An empirical Bayesian estimate of standard deviation for each measurement is used in calculating the standard error, which is a weighted average of the measurement standard deviation for each experimental condition and a global estimate of measurement standard deviation for each gene determined over thousands of arrays (Carlin, B. P. and T. A. Louis. 2000. "Bayes and empirical Bayes methods for data analysis," Chapman & Hall/CRC, Boca Raton; Gelman, A. 1995. "Bayesian data analysis," Chapman & Hall/CRC, Boca Raton). The standard error is used in a t-test to compute a p-value for the significance of each gene expression change. The coefficient of variation (CV) is defined as the ratio of the standard error to the average Log₁₀-ratio, as defined above.

Example 2

Preparation of a Chemogenomic Dataset for Late-Onset Renal Tubule Injury

[0115] This example describes methods used to prepare a chemogenomic dataset (i.e., a positive training set) for use deriving a signature for renal tubule injury (i.e., late-onset nephrotoxicity).

[0116] Overview

[0117] 28-day repeat dose studies were conducted on known nephrotoxicants. Doses were chosen that would not cause histological or clinical evidence of renal tubular injury after 5 days of dosing, but would cause histological evidence of tubular injury after 28 days of dosing. Animals were assigned to groups such that mean body weights were within 10% of the mean vehicle control group. Test compounds were administered either orally (10 ml of corn oil/kg body weight) or by intra-peritoneal injection (5 ml of saline/kg body weight). Animals were dosed once daily starting on day 0, and necropsied 24 hrs after the last dose following an overnight fast on day 5 (n=5) and day 28 (n=10). An equivalent number of time- and vehicle-matched control rats were treated concurrently. Likewise, a large set of short-term (day 5/7) treatments that would not cause renal tubular injury (i.e., negative control data) after sub-chronic dosing conditions were selected from the chemogenomic reference database in-vivo studies described in Example 1 (above), to complete the training set. This assertion of the absence of nephrotoxicity for these compounds was based on thorough evaluation of human clinical studies curated in Physicians Desk Reference (PDR) as well as peer-reviewed published literature. Lastly, these treatments did not cause histological evidence of renal tubular injury on day 5/7. Appropriate time and vehicle-matched controls for these negative treatments were also derived from the reference database in vivo studies described in Example 1.

[0118] Compound Selection and Dosing

[0119] To derive a signature predictive of renal tubular injury, it is necessary to first define both nephrotoxic and

non-nephrotoxic treatments from short-term studies devoid of tissue injury that can be used to model the early transcriptional effects that will be predictive of late-onset toxicity. To empirically confirm the late-onset nephrotoxicity of the positive treatments prior to inclusion in the training set, 28-day repeat dose studies were conducted on 15 known nephrotoxicants in adult male Sprague-Dawley rats according to the in vivo methods described in Example 1.

[0120] In addition, 49 short-term (day 5/7) compound treatments that would not cause renal tubular injury after sub-chronic dosing conditions were selected from chemogenomic reference database (DrugMatrixTM) to complete the training set. This assertion of the absence of nephrotoxicity for these compounds was based on thorough evaluation of human clinical studies curated in Physicians Desk Reference (PDR) as well as peer-reviewed published literature. These treatments were experimentally confirmed not to cause histological evidence of renal tubular injury at the time of expression analysis.

[0121] Doses were chosen that would not cause histological or clinical evidence of renal tubular injury after 5 days of dosing, but would cause histological evidence of tubular injury after 28 days of dosing. This time course of injury was significant to deriving a predictive signature since the presence of injury on day 5 would bias the signature towards a gene expression pattern that are indicative of the presence of a lesion, rather than identifying gene expression events that will predict the future occurrence of the lesion.

[0122] The compounds and their doses are listed in Table 2.

TABLE 2

64 in vivo compound tr	64 in vivo compound treatments used in the training set.											
Compound	Dose (mg/kg/d)	Time (d)	Vehicle	Route	Class							
4-NONYLPHENOL	200	5	Corn oil	PO	Nephrotoxic							
AMIKACIN	160	5	Saline	IP	Nephrotoxic							
CADMIUM CHLORIDE	2	5	Saline	IP	Nephrotoxic							
CARBOPLATIN	5	5	Saline	IP	Nephrotoxic							
CISPLATIN	0.5	5	Saline	IP	Nephrotoxic							
COBALT (II) CHLORIDE	10	5	Saline	IP	Nephrotoxic							
CYCLOSPORIN A	70	5	Corn oil		Nephrotoxic							
DAUNORUBICIN	4	5	Saline	IV	Nephrotoxic							
DOXORUBICIN	4	5	Saline	IV	Nephrotoxic							
GENTAMICIN	40	5	Saline	IP	Nephrotoxic							
IDARUBICIN	4	5	Saline	IV	Nephrotoxic							
LEAD (II) ACETATE	2	5	Saline	IP	Nephrotoxic							
NETILMICIN	40	5	Saline	IP	Nephrotoxic							
ROXARSONE	11	5	Corn oil	PO	Nephrotoxic							
TOBRAMYCIN	40	5	Saline	IP	Nephrotoxic							
6-METHOXY-2-NAPHTHYLACETIC ACID	360	5	Saline	PO	Non-nephrotoxic							
ACARBOSE	2000	5	Water	PO	Non-nephrotoxic							
AMPRENAVIR	600	5	CMC	PO	Non-nephrotoxic							
ANTIPYRINE	1500	5	CMC	PO	Non-nephrotoxic							
ASPIRIN	375	5	Corn oil	PO	Non-nephrotoxic							
ATORVASTATIN	300	5	Corn oil	PO	Non-nephrotoxic							
AZATHIOPRINE	54	5	Water	PO	Non-nephrotoxic							
BENAZEPRIL	1750	5	CMC	PO	Non-nephrotoxic							
BETAHISTINE	1500	5	Water	PO	Non-nephrotoxic							
BISPHENOL A	610	5	Corn oil	PO	Non-nephrotoxic							
BITHIONOL	333	5	Corn oil	PO	Non-nephrotoxic							
CANDESARTAN	1300	5	CMC	PO	Non-nephrotoxic							
CAPTOPRIL	1750	5	Water	PO	Non-nephrotoxic							
CELECOXIB	263	5	Corn oil	PO	Non-nephrotoxic							

TABLE 2-continued

64 in vivo compound t	reatments use	d in th	e training	set.	
Compound	Dose (mg/kg/d)	Time (d)	Vehicle	Route	Class
CLINDAMYCIN	161	5	Saline	IV	Non-nephrotoxic
CLOFIBRATE	500	7	Corn oil	PO	Non-nephrotoxic
CROMOLYN	1500	5	Water	PO	Non-nephrotoxic
DEXIBUPROFEN	239	5	CMC	PO	Non-nephrotoxic
ENROFLOXACIN	2000	5	CMC	PO	Non-nephrotoxic
ETHANOL	6000	7	Saline	PO	Non-nephrotoxic
EUCALYPTOL	930	5	Corn oil	PO	Non-nephrotoxic
FENOFIBRATE	215	5	Corn oil	PO	Non-nephrotoxic
FLUVASTATIN	94	5	Corn oil	PO	Non-nephrotoxic
GADOPENTETATE DIMEGLUMINE	125	5	Saline	IV	Non-nephrotoxic
GEMFIBROZIL	700	7	Corn oil	PO	Non-nephrotoxic
GLICLAZIDE	1500	5	CMC	PO	Non-nephrotoxic
GLYCINE	2000	5	CMC	PO	Non-nephrotoxic
INDINAVIR	1000	5	CMC	PO	Non-nephrotoxic
KETOPROFEN	20.4	5	Corn oil	PO	Non-nephrotoxic
LEFLUNOMIDE	60	5	Corn oil	PO	Non-nephrotoxic
LINCOMYCIN	1200	5	CMC	PO	Non-nephrotoxic
LISINOPRIL	2000	5	CMC	PO	Non-nephrotoxic
LOVASTATIN	1500	5	Corn oil	PO	Non-nephrotoxic
N,N-DIMETHYLFORMAMIDE	1400	5	Saline	PO	Non-nephrotoxic
N-NITROSODIETHYLAMINE	34	5	Saline	PO	Non-nephrotoxic
RAMIPRIL	1500	5	CMC	PO	Non-nephrotoxic
RAPAMYCIN	60	5	CMC	PO	Non-nephrotoxic
RIFABUTIN	1500	5	CMC	PO	Non-nephrotoxic
RIFAPENTINE	75	5	Corn oil	PO	Non-nephrotoxic
SULFADIMETHOXINE	1100	5	CMC	PO	Non-nephrotoxic
SULFAMETHOXAZOLE	1000	5	Water	PO	Non-nephrotoxic
SULFINPYRAZONE	269	5	CMC	PO	Non-nephrotoxic
TENIDAP	75	5	Corn oil	PO	Non-nephrotoxic
THIAMPHENICOL	1500	5	Water	PO	Non-nephrotoxic
TRANSPLATIN	0.5	5	Saline	IP	Non-nephrotoxic
VALACYCLOVIR	88	5	CMC	PO	Non-nephrotoxic
VALPROIC ACID	850	5	Water	PO	Non-nephrotoxic
ZILEUTON	450	5	Corn oil	PO	Non-nephrotoxic
ZOMEPIRAC	11	5	Saline	PO	Non-nephrotoxic

[0123] In Vivo Studies

[0124] Male Sprague-Dawley (Crl:CD® (SD)(IGS)BR) rats (Charles River Laboratories, Portage, Mich.), weight matched, 7 to 8 weeks of age, were housed individually in hanging, stainless steel, wire-bottom cages in a temperature (66-77° F.), light (12-hour dark/light cycle) and humidity (30-70%) controlled room. Water and Certified Rodent Diet #5002 (PMI Feeds, Inc, City, ST) were available ad libitum throughout the 5 day acclimatization period and during the 28 day treatment period. Housing and treatment of the animals were in accordance with regulations outlined in the USDA Animal Welfare Act (9 CFR Parts 1, 2 and 3).

[0125] Clinical and Post-Mortem Evaluation

[0126] All animals were monitored daily for clinical observations approximately 1 hr after dosing. For both the reference database studies (described in Example 1) and the sub-chronic study presented herein, gross necropsy observations and organ weights (liver, kidneys, heart, testes) were recorded for all animals following termination. Paired organs were weighed together. Body weights were recorded pre-test and daily thereafter for reference database (i.e., DrugMatrixTM) studies, and on days 0, 3, 5, 7, 14 and 28 for the sub-chronic studies. Terminal body weights were measured at necropsy and used to calculate relative organ weights and percent body weight gain relative to day 0.

[0127] Clinical Pathology

[0128] Blood samples were collected at necropsy from the orbital sinus or abdominal aorta under CO₂/O₂ anesthesia prior to terminal necropsy by exsanguinations and pneumothorax. A panel of clinical chemistry and hematology parameters were analyzed on a Hitachi-911 and a Baker 9000 instrument, respectively.

[0129] Histopathology

[0130] The right kidney was preserved in 10% buffered formalin for tissue fixation and subsequently embedded in paraffin, sectioned and stained with hematoxylin and eosin. Sections (5 µm thick) were examined under light microscope by Board Certified Pathologists for histopathological lesions. The left kidney was snap frozen in liquid nitrogen for subsequent RNA extraction.

[0131] Statistical Analysis of Animal Data

[0132] Treatment group means for body and organ weights, and clinical chemistry and hematology measurements were compared to the time-matched vehicle control group by Student's T-test. Significance was declared at p<0.05.

[0133] Microarray Expression Profiling

[0134] Gene expression profiling, data processing and quality control were performed as previously described in

Example 1. Briefly, kidney samples from 3 rats were chosen at random from each treatment and control group on day 5 for expression profile analysis on the Amersham CodeLinkTM RU1 Bioarray (Amersham Biosciences, Piscataway, N.J.). Log transformed signal data for all probes were array-wise normalized used Array Qualifier (Novation Biosciences, Palo Alto, Calif.), a proprietary non-linear centralization normalization procedure adapted for the CodeLink RU1 microarray platform. Expression logratios of base 10 are computed as the difference between the logs of the averaged normalized experimental signals and the averaged normalized time-matched vehicle control signals for each gene.

[0135] Results

[0136] A few treated animals showed histopathological evidence of early chronic renal nephropathy on day 5, including minimal to mild regeneration of tubular epithelium, interstitial inflammation, pelvic dilation, focal thickening of basement membrane and focal infarcts. Cisplatin induced a high incidence of mild tubular basophilia (4 of 5 rats), while both cisplatin and carboplatin induced a high incidence of karyomegaly (3 and 5 rats, respectively). Mild tubular dilation and proteinaceous casts were also observed in one lead acetate-treated rat. Although considered early signs of tubular injury, these mild and infrequent observations are unlikely to bias the signature since the large majority of the animals treated with the 15 nephrotoxicants were unaffected on day 5. Furthermore, the incidence and severity of findings indicative of tubular injury were markedly increased after 4 weeks of treatment relative to the day 5 time point.

[0137] After 4 weeks of dosing, all 15 nephrotoxicants showed evidence of degenerative changes of the renal tubules or early signs of tubular toxicity. Histological findings included tubular necrosis, dilation, vacuolation, basophilia, mineralization and cysts. These lesions were also accompanied by a higher incidence and increased severity of epithelial regeneration and interstitial inflammation, as well as granular and proteinaceous casts. A high incidence of karyomegaly was also noted for cisplatin, carboplatin, lead and cobalt. Consist with the tubular injury was the concurrent observation of hypercholesterolemia and hypoalbuminemia for a number of the nephrotoxic treatments. Although weaker than most other nephrotoxicants, 4-nonylphenol and roxarsone induced clear evidence of tubular injury on day 28. For example, proteinaceous casts, tubular cysts and mineralization were only observed in one roxarsone or 4-nonylphenol treated rat on day 28, yet these treatments did induce a much higher incidence and severity of tubular regeneration (4-6 rats) and interstitial inflammation (6 rats) suggestive of future tubular injury. Since the nephrotoxicity of 4-nonylphenol and roxarsone have previously been described (see, Chapin et al., "The effects of 4-nonylphenol in rats: a multigeneration reproduction study," Toxicological Science 52(1): 80-91 (1999); Latendresse et al., "Polycystic kidney disease induced in F(1) Sprague-Dawley rats fed para-nonylphenol in a soy-free, casein-containing diet," Toxicological Science 62(1): 140-7 (2001); Abdo et al., "Toxic responses in F344 rats and B6C3F1 mice given roxarsone in their diets for up to 13 weeks." Toxicology Letters 45(1): 55-66), and early signs of injury are apparent in the current study, these treatments were included in the positive class.

Example 3

Derivation of a Predictive Renal Tubule Injury Signature

[0138] Overview

[0139] The support vector machine algorithm was trained to classify experimentally confirmed nephrotoxicants from non-nephrotoxicants using the data acquired in Examples 1 and 2 above. A linear classifier (i.e., gene signature) was derived using kidney expression profiles from rats treated with 15 nephrotoxicants that induce renal tubular injury after 4 weeks of daily dosing, and 49 non-nephrotoxicants known not to induce renal tubular injury under subchronic dosing conditions.

[0140] Gene Signature Derivation

[0141] To derive the gene signature, a three-step process of data reduction, signature generation and cross-validation of the predictive signature was used. A total of 7478 gene probes from the total of 10,000 on the CodeLinkTM RU1 microarray were pre-selected based on having less than 5% missing values (e.g., invalid measurement or below signal threshold) in either the positive or negative class of the training set. Pre-selection of these genes increases the quality of the starting dataset but is not necessary in order to generate valid signatures according to the methods disclosed herein. These pre-selected genes are listed in Table 3.

TABLE 3

	7478 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #
NM_012939 NM_012657 NM_012848 U67914 AW915240 BF415939 L18948 NM_017250 AF150082 AI511090 AA859352 NM_017270 M63282 M35992 AB009636 X59132	AI180253 J02657 NM_012764 AB040031 AA818643 D38381 X83231 AB043981 NM_017288 U22520 BE113181 AB013732 D50671 AF202887 BE114586 AJ011607 NM_019126	AF139809 AF139809 AI7177121 D17310 NM_019308 X78997 AF055477 NM_013052 NM_019242 U75924 M96674 BE105381 NM_019322 AF034577 Z17239 AI029460 M11814 NM_013075	X63369 AI412259 AI011505 NM_012878 NM_019298 AB025431 M62832 AA849028 AA858817 AI175330 U16253 AW917537 AB042598 M81681 AI172112 AF306458	Accession # U27518 AF159103 D00753 AF290213 AI010583 AJ237852 AI410548 NM_013062 U56863 BF282409 U25137 D38101 AI407163 AW916143 NM_012698 AI575641 BF400833 J03863 Y13400
NM_013105 AF057564 BE109667 AF208288 NM_013068 NM_012682 NM_019233 NM_013197 AF151367 BF5555121 AI169311 NM_012738 NM_012786 BF522317 M26199 AB036792 AW143005 NM_012498	U04317 AJ276893 AI233740 BE100918 AF053312 AF044264 NM_012633 AB032419 NM_012810 J03734 J02635 AA997397 NM_012551 M22899 NM_017289 AF144756 M34052	AI171219 BF405468 NM_019348 AW920818 BF3995598 NM_019128 AI412261 X06827 AF199333 M74716 NM_017014 K03501	NM_012603 U66707 AI236696 BE109861 X05884 U94708 AF014503 J02643 AF058786 BE109018 NM_012803 AW916301 BE113155	NM_012639 AI236611 AF120275 NM_019286 AI009597 AW915049 NM_012567 AB000215 AF254802 AW141051 BF403190 NM_017123 AF227439 BE107840 U97146 AA893596 AJ001713 AI180010

TABLE 3-continued

TABLE 3-continued

	7/70	need to deale	DTI alarete			7470	need to deed	DTI alarete	σ.
	/4/8 genes	used to derive	RII signature	<u>s</u>		/4/8 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #
BF283270	AF112256	L19031	AW142962	NM_017215		NM_019339			AI102047
BF387347 AA891470	BE112719 NM_012735	NM_013086	AI409934	AI178784 BE112216	U31866 AI412108	BE108896 AF249673	J03093 NM 012588	AW918169	AW918050 L14323
NM_012881		AJ242926	NM 019344		BF285185	AI171162	Y00090	U66470	NM 017180
AA925167	AA901342	AI412418	L05435	NM_020087	BF556736	AW523849	AI228970	J03026	AW918529
NM_019295	X76723	AJ011035	NM_017279	AA800292	NM_012627	BF400832	NM_019326	AI136740	AW921215
AI234119	AF093567	M33936	NM_012614		AF295535	AA849743	AI454612	NM_017167	
NM_017354		X01976	AW143537	BE109691	NM_012825		BE107069	AI716512	AF148324
D87351 AF285078	NM_019310 AI233888	D89731	AI007992 AI008376	J02752 NM 012806	AI169596 AJ131563	NM_012842 U07971	AF15/016 AI411412	NM_013413 BE107234	AI576621
BF405086	NM_012879		AI012611	BF405917	M16235	AW251791	AI556066	BF550033	X53427
U61729	AI105410		NM_013217			NM_019204		BF563113	AW144705
BE105137	AA850034	NM_012870	U49066	AI010917	AW915996	U12309	X14159	NM_012851	AJ132008
NM_017259		AA819103	AF015304	NM_012533	BF283556	BE095878	AF198442	AA894092	AJ133104
BE113157	AI176677	NM_012757		BF401614	BF413176	BF282961	AW913932	BF283631	AW143091
AI574903 L17127	NM_012963 BF420018	AF063103 AF312687	AI137819 AW252871	D90109 BF542912	U41453 BF402407	NM_021691	NM_017074	M63122	BF556210 BF562701
AW914342	BF283381	BE111688	NM_012580		AF086630	BF405035	U33500	AI411995	U81186
AB012721	U57097	NM_012720		AI172281	AI407719	AA955213	AI045288	AW913986	AI232183
BF403552	BF416240	AI103158	AI603128	AW917780	NM_012938		AI101323	AW919092	BF411166
U80076	NM_012565		U15425	AW917985	BF398155	NM_017275		BF284803	M27223
U59245	AB005900	AA998157	Z17223	M15882	AA817877	U07560	AA817798	AA965057	AW917546
AI598399	AF111268	AW251703	AA946230	BF284124	AI172302	BE109271	AI230339	AF016297	NM_012835 BF394332
M94454 NM_021693	BE113285 BE113307	NM_012584 BE099881	U55995	AW915415 AW523614	BF562755 AF029107	NM_021746 M19651	AA848499	BF284475 NM 013057	NM_013176
AI176739	BF388223	AA848355	X87107	AI407487	AW862653	AF007212	AA892366	AW915287	X77797
U48596	BE098827	AF158186	AF068268	M84416	NM_012779		AI406538	BE113142	AA893184
AI412099	M58587	Y00065	U20796	AI180421	J02627	AW915613	U17967	AI168968	AF155196
U46118	U10188	AF133037	U41663	AW142880	X97477	BE119628	AI408557	AA875301	AF171936
AF027331	AI144646	AW920606	AW434178	BE113060	AF038591	BF285565	AI235942	AA964535	NM_012998
NM_012829 X15741	NM_017117	NM_017195	NM_017019	AI101117 BE282706	M94064 AI535126	NM_012908 AI170799	NM_012750	AW921399	U40064 AI764464
U44091	X94186	AI598316	NM_017208		AI059223	AW 144399	AW251848	AW914178	BF285034
AB017820	AF009329	AF109643	BF393825	X89603	AI234024	U12623	NM_012676		NM_012561
AF121670	BF284899	AI411981	AA800341	X68878	AI599016	AF009133	X65747	AI175457	AI012250
NM_013060		NM_019230		AI412460	NM_017113		AI406941	AI410352	AI408580
NM_013005 NM 012606		NM_017331 AI071251		NM_012833	AW917133 AA851926	Z78279 BF391604	AI236771 AF077354	BE120339 BF286916	AW144039
_		AW143506	AI555029 AI407201	AA945100 BF281697	AB020520	BF566679	BE109531	AA892778	AA800782 NM_019258
AI233903	AA892299	AI408713	AI411941	AA850910	AI102591	AF010293	BF551331	AW535307	BF283056
BE115621	AW921456	U26686	AF154114	U21871	D14015	D12769		NM_019289	
L27843	AW917933	AW915739	NM_021869		NM_012489		BE095859	AI176591	X70223
L29259	BF281701	NM_017097		AF089825		NM_017327		M26125	AA850347
Y18567 BF287903	U75402 AW915454	AW144649 J03886	NM_012618 AW917460	AII /1800 BF396132	AW915453 AA944169	AA818759 AJ222971	M18340 M90661	X73371 NM_017222	AI176836
NM_021836		AF184983	U67082	AI176814	AF013144	AW 252812	AB033771		NM_012845
AI111796	BF395192	BF414043	X84039	NM_013064	AF169636	BE116233	AF111181	U26033	Y09945
AW917212	AF105368	D83231	NM_012597		AI071412	BF406291	AI237075	AB007689	AW918231
AI010950	BF283340	AI227912	AA819832	AW914004	AI411400	D10699	AI454923	AF021854	NM_019309
NM_012771		AI408286	AI111954	NM_017115	D88586	AI104278	AW919284	AA800719	BF389519
NM_017011 X81395	NM_013008 U39943		AI716469	AW523875	NM_012964	X03369	D85760	AI172189	BF41148 M13646
NM_012794		BF288765 AA817759	BE105618 L19656	AW919125 M35297	BF407456 NM_017027		NM_012668 AI044740	NM_017089	
U44845	AW919210	D16237	NM_017261		AA850541	L19341	AW528864		AW251324
X83399	AI556458	AA892049	AF155910	AF036959	BF564217	NM_021747		AW524733	AB000776
NM_021763	AA925375	BE109730	J03627	AF041374	AA858862	AW917544	Y18965	AI409871	M80550
AI008409	BE106971	BE117330	L36459	AI137683	AW142828	U05675	M17412	BF288073	NM_012625
D50664	AF179679	NM_012621		AI412889	NM_017335		AF248548	NM_017310	
NM_017122 AI172222	BF396293	AW520812 AA800587	BE113132 NM_012940	AJ131848 AJ224120	NM_013106 U30290	X13016	AI235546 L25527	U39546 U61266	AA851370 AW916826
BF389915	Y17606	AF193014	NM_019358		AI010251	BF403184	NM_019904		NM_013226
AI549393	U35371	BF282980	AI008390	U09228	AI012235	AB000489	AI101924	NM_013043	
Z30584	U32679	AA858900	BE120346	AF001417	AW143771	AF136584	BF549650	AB023432	BF417292
AF102854	AW526005	AW915775	X12355	AF286595	BF549490	AI407141	AB039663	AW524724	D86383
AI237640	BE109016	AF277452	AW142290	AI177015	AF009330	AW434228	AI234852	BF412073	AI176773
L06821 X14788	NM_020084 AI408348	AI102884 AW919995	BF283610 X91234	AW862656 BF523561	AW525762 AW919683	BE106791 BF396602	AW921038 NM_012701	BF405050 AA859796	AW914097 AA799676
AW918179	M37828	AW919993 L46791	AJ1234 AJ137339	L15453	BE117335	BF524971	Z34264	BE349698	AF08797
AI716265	U15098	AF104362	BE108873	AA850740	U22893	AI170400	AW918222	BF284692	BE116886
BF551328	AI144797	BF415024	BE113252	AI179990	BE120016	AI411304	AW919395	X99723	AB027143
BF554744	AI176553	J03583	AI716560	AW918006	AA892250	BE103975	AB027155	AI230728	AI007666
D49977	U65656	M26744	NM_013139	AF172446	AF055286	NM_021595	AI406747	AI236770	AI406697

TABLE 3-continued

TABLE 3-continued

	7478 genes	used to derive	RTI signature	<u>s</u>		7478 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #
AW529723	BF558694	AW915643	AI406948	NM_012522	AI009609	BF285109	BF558592	BF282370	BF410020
AW913917 NM_013089	D87336 M15797	AW917481 BE349770	BF404362 M95058	M83209 AI176842	AW523642 AA901337	BF400636 BE117156	U17565 BF415786	NM_017065 NM_017276	
U66471	AW918776	BF550231	AA818471	L04527	AI145991	BE120810	NM_017147		NM 021585
AA955786	AI409218	U63740	AI176476	NM_019168	BE113179	NM_013074		AW253947	AA957047
AI012434	AA891839	X55995	AI406342	X53724	BF398114	AI072384	NM_019291		AF188608
AI575699	AF237778	AF327513	AF247451	AW525211	M38060	D16817	AF182949	AI227945	AI169105
NM_021587 AF118651	All 78875 AW 523888	AI177140 AB001075	BF553139 NM_016998	AA866351	U89608	NM_012679	AW141292 AI168941	NM_012678 BF418630	AW918/17 AA850896
AF118031 AI385364	U61373	BF556958	AW920082	AF297118	NM_012920 AF063102	AA943576	AI411510	NM_017337	
AI579216	AB043892	U93880	AW535233	AI406968	AF149118	AI409040	AA899704	NM 021842	
NM_017271	AF045564	AI171653	BE099774	BF404452	AI011501	AI598976	AA945761	BE110514	AW917461
BE100748	AI011757	BF399489	L05175	NM_017185	AI235960	U60282	AI406809	AF184921	NM_012923
BE106832	NM_017090		NM_017029		AW916210	AA945869	AA874859	AI176970	AI412936
U62667 AF082535	NM_017359 X85183	M55601	AF214733 X89968	BF393972 BF398716	AW917663 BE115280	AW918611 AF190458	AI102026 AI176993	BE115041 AA817769	BE098743 U94709
BE121120	AF013241	NM_012505		NM_012530	BF413244	NM_019153		AA955926	AF223677
U23407	AI007919		NM_013036		U15408	AA946474	U10357	AB012933	AW529756
AI409065	AI111579	AA799503	NM_019375	X78606	AA955579	AF007818	AA799664	AI172465	M35106
M77479	AW531805	AF072935	AW144670	D88666	AJ293617	X66539	AW913998	BE104266	NM_021775
X04644	AW921168	BF398063	BF282712	BE107169	L19927	C06844	X69716	BE116564	S57864
BE116867 AW917160	D85035 AW918417	AI235674 L13600	BE107427 NM_012942	AW916745 D13927	X99470 AA819481	AI070270 AW520754	AF118816 AI234012	BF562675 NM_013003	AA955206 AW532179
NM_017135		NM_017193		NM_021688	AI172146	AW918029	BE108973	AI599126	NM_019620
AA945724	AF080106	AF199504	BF410183	U42975	AW915254	NM_019369		U89873	AI172057
AW143008	AW434109	AI102073	AB016425	AF011789	AW920478	X95507	NM_019262		AW252815
BE105967	BE119862	AI105049	BF525016	AI170249	BF550822	AF022085	NM_019282		BF285046
BF551345	NM_016987		D88250	AW529672	U39571	BF284809	U92010	D83036	AA848826
BF554752 U06755	NM_017292 AW918564	BF542426	U52102 AI412423	BF283390 BF561659	U75920 AA893172	M34643 X76168	AF003835 AI111802	AI178527 AA943995	AI175508 U58466
U10303	NM_012959		BE095842	AB035306	AF314540	AI170766	AI111602 AI230699	AI406646	AB036421
AI410546	AI179101	AA900654	AW914408	AW921109		BF289272	D12678	AI412304	AI104545
L11002	AW920764	AI009727	BE103689	NM_012880	BF548454	BF389876	AA848311	AI599520	BF404935
AI169317	U64451	AA900261	AA819268	NM_021846	J05266	BF406752	AF286534	BF397840	AB017260
BF284190	AI502229	AI101181	AB011365	AW915558	AA800222	NM_012549		AA801208	AB021980
BF409783 Z14030	AB003587 AA800260	AI102732 BF395781	AF021348 AF192366	X3015 AA891535	AB011533 AI103924	NM_016995 U95368	AI234142 AI236084	AA849975 AI045904	AI012231 BE120595
AI237403	AB011528	BF563517	AW527151	AA892500	AW529808	X60370		AI411580	AW533822
NM_019192		AA818353	BE098709	AW524433	AW916093	X74815	AW917823	M55049	BE106398
X04310	J00705	AI176497	BF557244	AW916447	BF282239	AA800184	BF406312	X78461	AA817812
AF014827	U44750	AI227885	BF567710	BE107464	BF420074	AF306457	NM_012661		BE106693
AA943824	AI179459	AW254369 D25233	NM_012501 U89514		L29232 NM_017144	AW144456	NM_012967 AA866419	NM_012563	
AF186469 AW144673	BF412037 D13963	L14617	U91539	NM_012838 NM_019370	NM_017144 NM_019214		BE119802	AI413060 BF408325	BE108277 AA944061
BF290076	M60388	NM_012532		AI169655	NM_020096		AI104484	X06338	AA945771
NM_021868	NM_013134	NM_019283	AI234149	BE106663	X58465	AF025670	AI144644	AA892325	AB024930
U74586	X05883	AB024398	AI171288	BE109059	U42627	AW921149	AI411971	AF207605	AI010965
L35921	AA891690	AF151373	AI716500	NM_012715	AA925353	AA944556	AW435041	AF263368	AW143130
AF133731 AA997881	AF251305 AF163569	AI717736 BE113224	AW528865 BF288208	U73525 AA850319	AF059258 NM_021682	AI010281 AI011455	BF556879 NM_012693	AW916182 BE419280	AW523899 BE111805
AW916609	AI599801	BF420610	NM_012884		AA799526	BF282149	AB017696	BF556833	BF282313
BF282415	BF413513	AF230645		AI229630	AW25115	BF403319	BF282034	NM_019293	
D13121	M55636	AI178171	AW144339	AI716086	AW915925	BF416794	D85100	BE101619	AI010948
U44129	NM_012667		AW533508	BE106513	AI177397	D10926	AF100960	AW144663	AW251681
AF021343 BE097240	NM_019240	AW434520 NM_019186	AW917427	BF401626 BF404514	BE120725	AI104478	AJ011811 BE418860	BF412296	D70816
AI179988	AW916836 AI103918	U23377	BE101480	BF404514 BF415760	BF564940 AB002801	AW144344 AW534329	BF418869 NM_017235	BF416115 X69834	NM_017038 AA944518
D90404	AW531368	AA848451	BE107465	BF420628	AI009657	U59486	AA946430	AI235784	AI411149
NM_017047	BE113242	AF067727	BF549703	BF556874	M65148	AA799993	AI104125	AW535151	BF558676
U47280	BF285921	AA 998660	AF189019	NM_019143	NM_019216		AI144863	AW919320	Y14933
AI412625	NM_017340		BE108837	AW141939	AA892770	AI704771	AI409756	BE104321	AA850509
AI411021 BF416285	X70871 AI101396	AA818392 AI071243	AI010272 L22079	AW434239 AI030179	AI556534 AW914277	AW144347 AW913888	AW140637 AW918085	AB002111 Y18208	AI232565 AI407560
AI172196	AI101396 AI113104	AI071243 AI177050	AA894335	BE110739	BF248137	J05499	BF283003	BF420720	AI501497
BE116946	BF281931	AW435011	AF036548	NM_012553	BF563117	AB035201	BF408444	NM_017204	
BF398367	U14907	AW918443	AF157498	NM_012689	AW914809	BF404557	AI412054	AB052170	BF551283
NM_012999	AW253265	BE106598	AI105080	AA818796	BF287191	M95768	BF282646	AI598407	D49836
BE108670	M31363	BF411113	AF320509	AI237118	BF414010	X06564	BF396955	AJ132846	M18467
BF287788 BF404478	U30789 AF038388	AA965219 AI555002	AI012574 AW532606	AI639168 BF282476	BF566546 L22191	BE097587 BG153269	BF419234 BF556836	AW141993 L19118	NM_013091 AA891734
AW141130	AW526289	BE098212	AW915048	BF283760	L36884	AA891742	AI045635	NM_013148	
NM_019306		AF076183	BE108327	BF398053	L43592	AF030423	AW144346	NM_019155	

TABLE 3-continued

TABLE 3-continued

	7478 genes	used to derive	RTI signature	<u>s</u>		7478 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #
NM_017143		AW917390	AI105101	AW917475	NM_017212		AI599077	AA850490	X57523
U06436 AA850333	AJ242649	AW920271 BF281754	AI599133 BE107503	U50185 NM_017068	AA800476 AI231210	BF396180 BF399385	AW526320 BE116976	AW251878 AW531902	AF044201 AW914119
AI169058	AW441131	J03190	BF419138	NM_017068 NM_017260	AI231210 AI235446	D00403	BE117878	BF392577	AI102739
AW143818			NM 012869		AW531735	NM_017040		AI044865	AW920722
U85512	U69279	NM_013167		NM_017322	BE111795	AA850725	AF065161	AI059079	NM_013070
X07365	X60789	AB011529	AW535377	NM_017332	BE349699	AB041998	AF150091	AI454134	AA800062
AW525184	BF284889	AW531919	BE109179	U69702	BF284919	AW144170	AI231716	AI454913	AA818952
D12516	BF398564	AI228528	M37394	U77697	BF396319	AW920501	BE113635	AW919929	AI012608
AA800029	NM_013012		U13396	AF004218 AF058787	AF087696	BE108882	BF546202	D16479	AI137286
AI179538 AI227894	NM_017034 X82152	AI231808 AI412662	U17901 AI409316	D31873	AI009647 X74832	BF389244 BF408285	U04738 AI012120	NM_012545 NM_012555	
BE116816	AI232098	AW524460	AW917557	AF043345	AF002281	NM_019386		Y00480	AI411501
NM_012739		AW917674	BE113545	U05014	AI013919	Y17048	AW253043	AF068202	AW915713
X59037	AF097887	BF287209	BF405134	AA891447	AJ295748	AA944314	AA818438	AI406502	BE105713
AF043642	AI235923	NM_019257	AA819488	AI412169	AW144385	AI172618	AI103634	AW253742	BE107247
AI231761	BF405417	AB033830	D88450	U55816	AW252087	AI180353	AW918605	AW919881	U19967
BF559190	D13962	AI409182	AA866477	X55969	AW920527	AI235467	BE115947	U62326	U35775
AF259981	L13445	BF289240	AI102735	AA848951	BE111888	AI407999	BE116569	X57764	AF037071
BF396614	NM_019205	BF412673	AI227769	AI071288	AA875143	AI547463	BE118450	BE112952	AI009599
D88035 NM_012972	AF046886 AI598942	NM_017217	BE106459	AW530415 AW915057	AF202265 AI231444	NM_013137 U37058	BF397834	NM_013153 AA874924	AI172198 AI227748
AA799476	BE102426	NM_021856		BE107033	BF555116	AF157026	NM 012503		AW143395
AA891746	BE328941	AF111160	D25224	BF283353	D87839	AI009074	NM_012936		BE113449
AI170114	BF563262	AI104292	AA850987	BF404908	AW142811	AW526039	NM_019305		BF396079
AI407982	NM_019203	AI603127	BF400611	BF407675	AW915601	BE109164	X97376	AI178923	M17069
AW914850	X68199	X07320	D28754	L07315	BF288273	NM_012514	BE106816	NM_012974	U01914
BF550847	AF021935	BE113369	M80601	NM_017131	NM_012721		BE118122	NM_013215	
X66366	BF403098	BE120386	NM_012524		AA892791	AI411930	BF403852	U09583	AF315378
AI710683	AI011610	BF399607	X17037	AF106659	AA899304	AI602125	D63648	AI233769	BF284065
AW916287	NM_012841	NM 017325	AW526352	AF277900	AI176505	BF522885		AI385171	AI007922
AI235047	AF095576	U60096	AW915106	AI172003 AI411352	AI410099 AI709768	BF288063 BF565795	AA944006 AB017702	AI409259 BE108949	AI013500 AI170394
BF396316	AW918640	AF273025	AW914808	AW921986	BE111762	NM_012862		BE103916	AI556941
BF418597	BE119482	AI411225	BF283797	BE107674	BF282458	AA800455	BF415061	AI227890	AW141985
AA963234	U96921	BE107250	BF403853	BF285066	BF409042	M18028	BF389856	AI059493	AA800010
AF071003	AA818184	AW915650	BF555793	D49955	U06864	NM_012733	BF393807	AA942765	AA891944
AF182717	AI168986	BF285467		NM_012497	U08136	X83579	D16348	AI137297	AF000577
AI574743	AW915209	BF398403	U60835	AI233257	U50707	AA997412	NM_017316		AF017393
AI599339	BF286131	M29853	AI169629	NM_017075	AF245172	AF151982	U10279	AW252664	AF065147
AW914966 AW915217	U16655	NM_012725 X51991	AW140983 AW916321	NM_021740 BF413765	AF272892 AI103040	AI137817 AI411605	AF290194 AI170797	BF397998 AF021923	AI100769 AI170067
AW917653	AW919239	X57970	D50568	BF550866	AI103040 AI227907	AW143854	AI170797 AI179993	AI230762	AI170405
BE113656	BE109746	X92069	NM_012903		AI412317	AW 507078	AI713159	AW141128	AI236618
BF284918	BF558086	AA892522	AI228236	M21208	BF404556	AW525122	AW917818	AW534159	AW918017
NM_012499		AF160978	AW143256	AF036537	NM_013049		J00741	AW915055	BE096387
NM_013145		AI711516	BE107147	AI600221	U67884	AW918208	AW251204	AW915084	BF288129
BE107259	BF398684	BE119676	AF030091	AW915518	AW914944	BE109208	BE095833	AW916127	BF403009
BF283261	AI072493	AA996836	BF287099	AW917510	AW916786	BE113233	BE110525	BF404842	BF555971
BF400209	AF099093	NM_017124		BE097210	NM_012828		BF400666 NM 012763	BF418582	NM_019273 X02904
BF388763 BF550217	BF524872 NM_013001	AF228917	NM_021703 AA944542	BF420172	AA900400 AF306394	BF396948 BF549525	U67080	AA943742	AA800273
M55045	U61157	AW915404	AF001896	BF557300	AW533663	NM_019363		AI168952	AI102064
NM_012704		AW917389	AI010954	BF559836	AW915554	AF239674	AI009650	AI556408	AI171802
U73458	M83676	BE114160	AI012498	BF567496	AW917516	AW915161	AI555351	D50559	AI230430
BF557672	AF163477	BF283798	AI176695	D25290	BE110577	BE105872	AI711114	AA945320	AW917132
X57228	AI175555	BF567631	AI180420	M88709	BF281848	NM_012819		AF247452	BE108178
BF291167	BF416236	L02530	AI406310	U32314	D86711	NM_019237		AI070113	BE108857
AI599349	U20195	AF080468	AW531361	X16359	J03819	AI406821	BF399083	AW915174	BF397872
AW520770 DE555127	U31668 AA850037	AI176944	AW916592 DE112033	AA801139	NM_017077 AW142947	AF084576 AI060205	BF419406 NM_013115	BE111769 BE116370	BF543356
BF555127 D12498	AB020757	AI179372 AJ003065	BE112933 BF290834	AI407016 AI412627	AW 142947 AW 434045	AI000203 AI179609	NM_021744		BF543478 NM_012655
J05132	AI009371	AW143887	BF412769	AI555466	AI232337	AI1/9009 AI408442	AB043870	BF554877	NM_012033
U48249	AI231799	AW916474	BF555129	AJ011608	AJ277747	AW915550	AI598442	BF556614	U16858
AW918418	BF396682	AW917766	BF557296	NM_017181	AW918255	BE113312	AW253880	U29174	AA851728
BF281400	BF417187	NM_013071		AB004329	AW919873	BF387255	AW917588	X85184	AI137188
D16308	M63574	NM_017345		AF205438	BF286478	BF394261	BE112998	BF566580	AI177431
AW917550	AA799741	X71068	U17604	AI407017	BF388422	AA964824	U41853	AI176632	AI555341
BF400606	AA799751	AI411997	AA819306	AW143149	NM_012562		AA817841	AI178935	AI600037
BF547620 BF563077	AF106860 BF281149	AA800507 AA875261	AF005099 AW918548	AW918238 BF386716	NM_019211 U42209	AU231792 AW521352	AA849966 AF220760	AI406531 AW144006	AW526346 BE108174
M61726	AB019791	AA875261 AI231776	AW918548 X74549	U53855	AI412190	AW 521352 AW 917064	AF220760 AI102512	NM_017141	
14101/20	211017171	111231110	411 マンイン	000000	2 TT T T T T T T T T T T T T T T T T T	2 111 / 1 / UUT	Z 1102012	7.11.1	202121760

TABLE 3-continued

TABLE 3-continued

	7478 genes	used to derive	RTI signature	<u>s</u>		7478 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #
AI716642	AW918585	AI231088	X04240	BF282876	BF544951	AB014089	BF282271	AF090113	BE109596
BF555161 NM_012527	BE101088	AI598315 AI713206	AB023634 AI175008	BF557821 D31838	NM_017005 U77933	AII /5820 AI237593	D78482 AA893708	AF127390 BE113175	NM_017189 AA899898
	NM 020106		AI237657	L27081	AA998160	AI409747	AB001982	BE113372	AI410203
AB021645	X71071	BE117902	AI717113	NM_019167	AI176825	AW253010	AI716516	BF288088	AI705687
AI169116	AA955396	BF282238	AW252879	AA849738	AW535229	AW143285	AI172159	AI231206	BF409759
AW919190	AI176056	BF410755	AW532489	AI235219	AW915543	AW916783	AI409738	BE108849	AA859585
BE098463	AW916119	M22323	AW916092	BE107395	BF284679	AW917522	AW915402	BF389882	AF109393
NM_012590		M81784	BF564219	BE108776	BF285207	NM_012685		BF550292	AI009274
NM_019364		Z48444	NM_013065 U49057		BF392605	NM_012728 X53003		AI010241	AI013361
U11038 AA850505	AI179901 AI407827	AA799832 AF096835	AA800483	BF550270 AA859010	BF398046 AI008125	AF098301	AI556488 BE115058	BF558976 AW915795	AI013475 AW525285
AA892818	AI716077	AI008701	AB042407	AI600035	AI172184	AF199411	BF283742	AA894080	BE103518
AI177408	BE098955	AI176212	AI170313	AI716255	BF550875	AI231432	AI103682	BE097615	BE114137
AI227612	BF388440	AI176625	AW252251	BE102621	AF009603	AI236772	AW534533	BE108899	BF289044
AI412150	AF277903	AI180275	AW917256	BF284716	AI009591	AI408517	AW535909	BE113057	J05030
AW523647	AI058960	AI412673	AW919062	BF397805	BF392344	BE126739	BE098266	BF396082	U77038
AW917504	AI409077	BE113146	BF389884	BF400811	BF404539	BF288138	D90166	M35052	U95727
BE095865	AW141921	BF403136	BF396282	L18889	AF000423	BF396678	U07201	M84009	AA943100
BE103444	BF285451		NM_012671		AI010267	BF558506	AA891834	NM_017178	
U68168 AF053317	AI229849 AW918000	BF283898	BF401593	NM_013223 U40819	AI101199	BE111625	AA997458	NM_019379	
AI169878	BE120513	BF284303 BF555924	BF413556 NM_012916		AI231787 AI715452	BF558467 NM_017152	AI044229	NM_021576 AF267197	AI143339 AI172450
AI109878 AI178796	BF419158	NM_012795		AA850358	AJ002940	BF419635	AW921797	AF276940	AI175031
BE107485		NM_013135		AI011497	AW525049	D21799	BF412792	AI454081	AI234810
AA943552	NM_017365		AI230723	AI070137	AW919666	AI069912	D13127	AW141938	AI408705
AA943564	NM_019219	U43175	AI598405	AW532074	BE107334	AW914758	D89514	AW918816	BE099401
AF017437	NM_021771	AJ000696	AW143111	NM_017186	NM_012636		U55192	BE103359	BE120608
AI172175	U30381	BF548957	AW434242	NM_021759	NM_019284		AI045819	BE118465	BF550426
AW142913	AA892370	BF549697	AW920179	U50194	AA997435	BF404603	AW144075	NM_017159	
AW143093	AB022014	L38615	BF406407	AA866426	AI236090	BF555890	AA945706	NM_017311	
AW528898	AF020618	AF226993	BF413396	AI408960	AI575940	M64301	AA945734	AI103954	NM_013103
AW915175 BF287814	AF059311 AF090867	BE112913 D88364	D00252 AA799400	AW142588 BE113966	BE349755 BF281802	Y00102 AA858509	AF106945 AF142629	AA819871 AF083418	AI103456 BF284887
U61261	AI172386	NM_013114		BE117883	U56732	AI105215	AF176784	AW918470	BF409560
AB009463	BF284171	NM_019255		AF075382	AF292116	AI237580	AI102248	BF551138	AI235238
AB032164	NM_017220		BE116768	AF087454	AW918457	AJ293697	BE095605	AA800701	BE109510
AF031483	AA925490	U49055	M77362	AI175048	BF408873	BF565344	BE121438	AF052042	BF525211
D29969	AW142307	AA998662	AA850728	AI407222	BF409812	NM_012581	BE329061	AI013104	AI172460
M34083	AW435429	AI104846	AI178761	AI599104	AA800241	AI170786	BF550271	AI407821	AI233875
AA848470	AB012233	BF406522	AW916628	AW141280	AF050159	AI231438	L31840	AI598402	AW916561
AF058791	AI175440	BF413631	BE098366	BF419602	AF313411	BF408022	X64411	AI599376	BE108405
AI176665 AI178491	AI409380 AW520758	D38072 U75916	BF281544 BF523059	L02896 NM_017262	AI317840 AI412209	AI598462 BF281386	AA998893 AI101490	BF285247 BF285980	BF282009 BF555349
AI178491 AI232898	AW 520758 AW 524559	AI227843	D16302	AI406693	AW919129	AI112622	AW915318	U68726	BF556162
AI233288	AW534383	AI411425		NM_021762	M20406	AI172033	AW915609	X78604	BF562149
AW532652	AW915350	AW531412	AI230918	AA955175	NM_021857		BF407740	X90710	NM_017241
AW915437	BE113217	BE115635	BE110626	AW915491	AF131294	AI409049	D00680	AI179119	U26397
U06099	BE116973	AI232321	AF169409	K03250	AI234849	AW533060	AF010131	AI411742	AA900983
U83112	BF282030	AI236624	AW535358	AA891859	BF415080	AW919094	X79860	AW142808	AA965117
AA965063	BF284994	NM_012984		AF106657	X58375	AW920600	AA943981	AA817907	AI171654
AB009999	AI009603	AI009623	BF282629	AI105272	AI102037	BE101138	AW916468 NM_017101	AI179443	AI177089
AF023657 AF135059	AI011034 AI011713	AI010235 AI179979	BF556943 BF563933	AI170757 AI233199	AW251849 AW527217	BE107520 NM_012674		AI578861	AI408686 M97754
AW144034	AI171480	AI599143	NM_019362		NM_017066		AF314960	NM_017213	
AW251238	AI232365	AJ245707	U55765	AI410700	AF007549	AF037199	AI008988	U78889	BE107747
AW254429	AI408357	AJ306292	AI103327	AI598467	AI008386	AI145784	AI233241	AA891790	AB006461
AW915944	AI716103	AW917197	AW144790	NM_017199	AI104546	AI177867	AW143117	AA925922	AI234008
BF401313	BE111634	BE103937	BE108865	U09793	AI176039	BF283802	BE101096	BF408391	AA944483
J05122	BE112615	BF545951	BF395678	AA851327	AI235512	BF396424	BE108272	BF525153	AF322224
NM_019135		BF556845	M11563	AA945202	AI407464	BF405032	L34821	AI407903	AI763565
U12187	BF408081	AA944161	NM_017194		AI549323	BF416249	AI177887	AW914881	AW916701
X03475	NM_012601		NM_017267		BE103152	BF283736	AI410438	AI407483 AI535483	BF282212
AA849719 AF268030	NM_017342 NM_021764		NM_017344 U71293	AI102063 AI232205	BE108583 NM_017099	NM_017021 AI138061	AI410822	AW433595	BF401710 J02997
AI009594	U75689	AI579023	AI716115	AW916594	AA817863	AI138001 AI412244	BE099629	BE108976	AA848338
AI136848	X60822	AW141186	AW141000	AW921320	AB030644	AW915966	AA801434	M59742	AI454466
AW143233	AF063939	AW916805	BE109756	BE114154	AB042887	BE105397	AA819679		AI555844
AW144502	AI179974	BF291161	BF284075	AI411071	AI103943	BF417391	AF084241	BE113624	BE098873
AW251339	BF398016	D14908	AI716289	AW142171	AI170377	U18942	AW915444	BF406637	BF395080
AW435110	BF525193	AI176016	AA956784	AW434007	AI179991	U75973	AW918431	U92803	BF414124
AW916721	D10665	BE109747	AF058714	AW915587	AW435010	X62952	BE098309	AA850785	BF546361

TABLE 3-continued

TABLE 3-continued

	11.	BEE 5 CON	imaca			111	DEE 5 con	imaca	
	7478 genes	used to derive	RTI signature	<u>s</u>		7478 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #
AW526079	AA848834	BF407209	AB020759	AW144002	BE101311	AF245040	L32591	BE100014	AW918999
BE108249	AA894259	BF407452	AF036344	BF419074	BE109604	AI136871	AI010234	BE109057	D90036
BE109637	AF022952	BF550795	AI137471	BF552916	BF289328	AI177706	AI233766	BE119961	NM_021684
BE113111	AI408852	BF555867	AI145625	M64711	BF393085	AI180454	AI716240	BF397933	AA800597
BF398605	AJ005113	D13871	AI172211	AA800210	BF551339	AI231601	AW254017	AF034214	AA892281
NM_012836		NM_019220		AA850498	L37293	BE108326	AW919336	AF190798	AI169225
NM_013216	BF549710	U72994	BE112948	AA893230	AI010342	BE115880	BF415023	AI010660	AI234095
BF389352	NM_012839	AI170827	BF283612	AF115282	AA851945	BF394214	J05029	AI170570	AI411077
AI071698	NM_021653	BE113005	BF284840	AI169619	AA943868	BF399328	NM_019385	AW526160	AI639139
AI175474	NM_021865		BF414261	AW531530	AA963282	L10072	AA799515	AW531675	AW433942
AW917280	AW536019	BF389478	BF522056	AW919429	AJ293948	NM_012592		BE111118	AA892483
BF551315	BE099732	BF412016	D13061	BE102505	AW143480	NM_012793		BE118222	AI104485
AB016532	D12770	U57362	U46034	BF419925	AW915268	AB042599	AI102877	BF555119	AI407945
AI230220	AW143273	AI012356	AA799661	M76591	BE107438	AF156981	AI176623		AI409108
AW915159	AW523874	AI169243	AA875055	NM_013063	BF556273	AI176323	BE098468	AA892339	BE095970
BE108853	NM_019180		AA943094	AF072124	BF559875	AI317817	BF550402	AB002406	BE101099
AA891830 AI411897	AA874838 AF228049	BF287768 BF396114	AF037350 AF244349	AI177645 AW918369	M23984 NM_012997	AI599641	NM_021770 AA819398	AF203906 AI010233	BE107434 D32207
BE110722	AI412591	U68544	AI180400	BE120038	AA892298	BE097245	AA946128	AI010233 AI175028	NM_013034
BE112999	AW434329	AA800232	AI603627	BF284819	AI029960	BE109513	AF151377	AI406667	X89963
D26179	AW914982	AI104857	BE095490	BF406693	AI409930	AA801206	AI177663	AI407482	AI231190
L06238	AW917734	AI105461	BE109529	NM_012578	AI716131	AF231010	AI412090	AJ242554	AI412736
NM_017050		AI230228	BE113119	NM 017353	AW526697	AI413033	AI412292	AW434419	AW433944
U03708	BF386111	AI412612	NM_020089	AI101475	BE100193	AW143939	BE102889	AW521367	AW917545
AW915834	BF397542	AW140530	AW140531	AI176781	BE108131	AW531093	BF408844	BF283772	BF283384
BF284693	BF549877	BF555370	AI176792	AI411194	BE113228	BF282636	BF564899	BF388772	BF420685
AA944036	AI172191	AI170769	AI236760	AI705731	BF567904	U48247	NM_012634		AA944568
AI102429	AI232217	AI170280	AI598324	AW141990	M81766	AA849715	AA924526	BF550302	AI072892
AI171775	AW528823	AI179677	BE107173	AW253902	Y08981	AF020046	AA944278		AI105210
AI406506	BF285991	AI410505	L46865	AW524517	AW144637	AI412580	AF094821	AB026057	AI236773
AW531891	BF565628 AI235353	BF403332	NM_012987		AI009167	AI600237	AW144383	AI172214	AI406363
BE107157 BF404868	AJ300162	AW142852 BF286955	NM_017175 AA817895	L26450	AI408865 AI412011	AW915560 AI598414	NM_017251 AA892364	AI235950 AW916305	AI408954 AW915466
D12771	AW918833	NM_012515		M34384	AW915292	AW915580	AB020022	BE100201	BF417386
AA893610	D90102	AI013928	AI010432	NM_020306	AW915499	BF550554	AF051561	BE105305	BF551118
AB038387	U87305	AI176626	AI169228	U15211	BE113053	D63665	AF177478	D17447	D14013
AI170859	AA892330	AI233205	AW534781	AA850551	BF282223	AW254190	AF323615	L02121	NM_012947
AI234035	AI407409	AW142713	BF390657	AF051895	AI169291	BF555084	AI071688	M20133	AA998435
BE105286	AW144331	AW142877	NM_017207		BF399098	NM_021754		M34253	AF080568
BE111776	AW915847	AW915294	AF090347	AI412323	AA946357	AI009029	AW915256	AW919017	AI045590
BF281438	BF557668	BF392695	AF030377	U65007	AI008952	AI227700	AI178647	AA875129	AI070591
BF404419	AA848342	BF397773	AI102519	X66842	AI103937	AI409145	AF072509	AA900046	AW915160
L36388 X86789	AA942695 AA955630	M32061 X62528	AI177143 AI232354	AB026288 AI717447	AI227742 AI411999	AW525288 BF396534	AI172156 AI176848	AA946441 AF139830	BF285089 BF547641
AA849782	AF020045	AA849497	AW522044	AW142440	BE109075	BF404409	AI170646 AI407459	AF205604	U93197
AA874906	AI137298	AB026291	AW917726	AW 527204	BF399614	J03753	AI411005	AW252550	AA924945
AI169368	AI179370	AI317813	BE106275	AW915676	M29295	NM_013006		AW916799	AF000942
BE109266	AI575703	AW918480	U53475	AI409024	NM 012665		AW252152	BE111887	AW535349
BF390003	AW141463	BE111685	AA850801	AW253750	Z83868	AA819812	AW916013	AI102290	BF558075
AI103914	AW143992	BF281285	AF012714	AW535136	AI101393	AI169599	AW916792	AI233162	AI411332
AI170783	AW918108	BF396317	AF146738	AW917211	AI547421	AI227919	BF387153	AA799789	BF285720
AI713210	BE105452	BF548520	BE103926	NM_017356	AW143757	AW919050	NM_019341		BF557889
BE098845	U19485	D50580	BE109586	AF061266	AW525128	NM_013076		AI102236	AA944526
BE102816	AA946490	U34841	BF396467	AI012352	BE108832	U09229	AA892271	AI411240	AB049189
BF283510	AB040807	Y17319	NM_019381		BF403323	AA945103	AI008961	AA799301	AI101322
BF391673	AF039033	AW918273	U72660	AI412018	BF407165	AB018546	AW918092	AI236816	AI102495
X56541	AF092207	BE121325 NM_012980	U83897	AI600031	BF555033	AF182946	BF282185 BF395777	AI409186	AJ277881
AA800172 AF327562	AI072958 AI178489	NM_012980		AW433866 AW917796	M58716 NM_017188	AI180081	BF393777	AI012573 AI172116	BF409313 AA818820
AJ238717	AW434972	NM_020976		BE107459	AB047002	AI410886	BF420629	BF282323	AI102873
AW918775	AA848526	X96488	AI145039	BF399791	AI232269	AW915104	BF557739	BF283075	AI179142
BE104931	AF063447	AF272662	AW143197	L31884	AW918541	BF407878	J03637	M69056	AI230778
BE119692	AF218575	AW144391	AW918637	AI012263	BF523077	BF414947	Y12009	AI105441	BF285078
BF283247	AI170251	BE099953	BF284879	AI233726	L11319	D10655	AI175375	AI407500	NM_012659
BF555980	AI235480	BF282288	BF565705	AI408104	M23601	M92042	AI230185	AI170752	U18650
BF564461	AW144226	BF282645	U11685	AI555237		NM_017231		AI172417	AI013775
NM_013058		BF413969	U13253	BF285079	AI411520	NM_019335		AI412239	AI411964
U48246	AI013913	AA874952	AB017793	BF417363	BE108919	NM_020073		D85580	BE109603
AI013699	AI137301	AW915060	AI230988	M84488	BF558507	NM_021847		J03933	BE114159
AI409741	AJ001184	BE104111	AI385140	NM_021997	NM_013090		BF419380	L27513	M86870
AW142955 BE096047	AW917946 BE112415	BF283001 BF284914	AI407991 AW434026	AA858786 AA894084	Z83035 AB002466	AA850736 AI407932	AI409032 AW144517	NM_012911 AF281018	AB005549 AI231193
DE020047	D1112413	DI 207714	2 X 11 T J T U Z U	2 № 102 TO OT	2 ED002400	2 MTO 1732	2 x ** 1 7 1 3 1 /	211 201010	111231173

TABLE 3-continued

TABLE 3-continued

	7.470	14 1 1	DTI :			7470	14 1 1	D.T.I.	
	7478 genes	used to derive	RTI signature	<u>s</u>		7478 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #
AI230056	AJ005425	AW525342	AI013788	AI385277	BF400873	AW142549	U84038	M22631	AI012456
AI410833 AI555566	AW143263 AW917908	AW914215 BE103434	BF399587 M81687	AI409841 AW915241	BF551369 AI176483	BE098806 AW144499	AI071470 BF555429	NM_012609 BF564158	AI137208 A1598988
AI553500 AI598648	BE106888	BF389721	AI413058	BF398378	NM 013042		AI598321	AA800290	BE112781
AI716218	BE111752	BF397663	AF069525	J05214	AA818571	AI706767	BE111696	AW434213	BF393577
BE101448	BF282437	BF411381	AI060118	NM_012818	AA943149	AW915737	L20822	AI231846	BF414252
BE102671	BF290638	NM_012846		AF019109	AI169160	AW918850	U08141	AI408197	BF558120
BE118605	AF016049	NM_017216		AI104376	AI411217	BF283053	AA800519	AW525033	Y17325
BF555974	AA892897	AW435310	AA875011	AI228233	BF282194	AA799614	AF016047	BF284076	AI105265
AA817722	AB015433	AW917572	AA891774	AI639162	BF401587	AB032899	AI233267	M36074	AI112074
AI233194	AI234830	BE108192	AA892554	AW917587	NM_021594	AI406853	AW527592	U60063	BE099063
AI408375	AW141787	U76997	AI715257	BE100208	AA891221	AW527606	AI071703	AI169278	BE101628
BE109600	AW143141	AA892567	BE113288	BE108905	BF556691	BE112252	AI145019	AA801230	BF549638
BF567692	BF548116	AA999042	BF551361	NM_019280	AB010467	BF416533	AI412626	AA892319	M97380
	NM_012571		AA892346	NM_019622	AI102685	NM_017246		AF065438	U40628
AF085693	AI175536	AI599031	AI234858	AA944162	AI177409	AA946356	BF558459	AI102139	AW914919
AI171807 AW528847	BF554895 BF563786	AW915803 BF396191	AI602172 AI172579	AI137972 AF192757	AI229166 AW918105	AB017711 AI178752	X16481 AI410901	AI236798 AW916666	BE107373 NM_017274
AW920802	J00696	M64300	AI717053	AI170933	BE113010	AI178732 AI599125	AW915787	BE113034	AB008161
L35767		NM_017187		AW 529588	BF281834	AW144760	BE108235	BF284695	BF288270
AW143336	AI236376	AW915800	BF283743	AW 530272	BF386665	BE108884	BE108381	U36786	BF397445
AW144084	AI407946	BF282620	AI144583	AW918408	BF394140	BF284699	BF549121	AF036760	BF416387
AW252169	AW144223	BF401275	BE102535	BF283302	X13549	AA956764	BF559056	AI177061	NM_013111
AW528454	BF283130	AW917258	AA849729	L20821	U82626	AI112512	AF051155	NM_019372	
AW915763	L19658	AI233133	NM_017200	U04319	AA943793	BE107281	AI412014	AA892300	AF121893
BF419241	AA801116	AI408930	AW913858	BF410753	AI105167	AI176713	BE107155	AF032872	Y15748
BF557396	AI011704	AW918153	AI012474	AF286006	AW144315	AI178763	BE109130	AI103962	AA800044
M58340	AW144504	BE109152	AI412614	AW252511	AI236054	AJ299016	Y08172	AI176002	AW144441
NM_013190		BF288288	BE109201	AI170384	BF389493	BF406240	AI598359	AW916151	BF414262
NM_021750		BF562779	BE109644	AI410837	BF400662	AF184920	NM_017151	BE111638	AA848367
AI175586	BF407799	U54632	BF281325	AA894189	AI230729	AI072236	AW528874	AB031014	AI178206
AI411060	U07181	AW915140	BF523098	AF119667	BE115551	AW917568	BF550453	AA892294	AI229655
AJ001529 BF282544	AB020504 BF283685	BE109575 AA899489	D30795 L09653	AF228307 AI234719	AI012951 AW917662	BE112921 M73714	AA893193 AF181992	AI177845 AI411497	AI406371 BE117002
BF408448	BF405110	AI111840	NM 017105		BE113247	AI102744	AW144745	AW253895	BF282296
AA851280	U58858	AI412967	BE103894	AJ001044	BE099563	AI232494	AW252105	AW915264	AW433959
AA944380	AI172274	AI575671	AA799981	BE107298	BF548170	AI233702	AW526756	AW916138	AF029310
AI176442	AF159626	BE100155	AA943811	AW916684	AA849756	BF284127	BF399124	AB020879	AI103375
AI237621	BE115860	NM_013055		BF389719	AI229596	BF405996	U87627	AI171276	AI176541
AI409180	NM_013028	NM_019246	AI236778	BE108876	AF158379	BF522695	AI104348	AI1712840	AI227815
AI410943	AW140640	AI231333	AW143201	AI411399	AI170263	AI412601	AI231785	AJ000347	AI411985
AI411979	BF393884	AW523114	AW254246	BE118972	AI234844	BF412389	AI411141	AJ292524	AW142847
BE108923	NM_013185		AW916618	H35082	AI639157	BF414338	BE110537	AW142931	BE113048
BF386302	NM_017024		X04959	L34039	AW915774	AA799576	M11185	AW144646	AI179335
D30035	AW916148	U06713	AA800199	AW143157	AI232784	AF296131	AF029690	BF403923	AA801136
NM_012586		AI105154	AA819716	AW533321	AW916344	AI385216	AI010722	BF420067	AA817945
X56228	BF285301	BE109143	AA946074	BF412594	BF408552	BE110949	AW252820	NM_019275	
AF003944 AI013474	D90035 AF220455	BF567763	BF396729	BF567585	AI233916	BF284939	AW914860	BE107208	AA850909
AI1015474 AI101500	AF220433 AI104326	NM_019208 U04933	AI104251	AW142367 BE121429	AI409258 BE098359	BF555949 BF564549	BF405883 L16532	AI103616 AW144313	AA891818 AI104296
BF284242	AW140537	BE096021	AI231564	BF407916	BF418913	L20900	AI012438	AW 529753	AI231812
U58857	M94548	BE113323	AI231789	M86235	J04112	Z16415	BE329046	AW915952	AW252855
AI029291	AA924352	BE121314	AW253339	AI009759	AA945604	AI229684	AI136513	AW918376	BE103222
AI170751	AW916619	BF407511	AW524478	AI407545	AB017544	AI406527	AI169330	BF404589	BF288776
BE112253	AW917712	NM_017079		AW918385	AI170948	AI409951	AI171772	BF410846	BF394038
NM_021848	BE108877	NM_017174	BE117114	BE101157	AW143214	BE098713	AI407001	BF419489	BF397229
AI071187	BF284713	AA849031	BF404464	AA799507	BF283454	M31788	AI548694	BF567996	BF558902
BF405880	AA799636	AA859343	BF563403	AA818132	BF523555	U14533	AW920624	X62322	NM_012804
BF548241	AI407904	AA943765	D50564	AI102046	U70825	AI178912	BE115875	AI044638	NM_016988
M93271	AW254590	AI175728	NM_017033		AI176331	AI408244	BF397523	AW529960	X94351
NM_019222		AI228548	AI101900	AI172271	AI013800	AI704755	BF408216	BE095474	AA893217
U82623	BE104941	AI230073	AI413051	AI230110	AI412560	BF282119	BF558463	BE108346	AA943578
AI410415	BF393950	AW433847	AW917849	BE102814	AW914984	BF392959		BE110542	AB028934
AW142953	X87885	AW915824	BE100016	BE118552	AW919694	BF409371	AA849991	BE110542	AI172177
AW434978	AI169383	BE098021 U79661	BF404932	BF404472 NM 013221	BE113234	NM_013166		NM_019299	
BE100035 BE108780	AI412413 NM_012528		BF416377 U23443	NM_013221 NM_021592	BE113330 BF398543	X87106 AI013041	AA894233 AI010295	AA893811	AW919497 BE111972
	AI412230	AII /5 /62 AW918595	O23443 AI408984	AA944463	M57299	AI013041 AI172285	AI010295 AI011448	AI178257 AI711105	BE111972 BE118440
AA849752	AW525071	BF281282	AI408984 AI411771	BF281215	NM_016986		AI011448 AI229529	AW142280	BF283418
AB003042	AF061947	X78949	AF065387	AA894318	NM_017153		BE105699	AW 142280 AW 915107	BF420144
AI236861	BF389157	AA963096	AI176933	AI009656	BE101101	AW915155	BF413204	AW915928	BF550580
BE099603	AI008969	AA998971	BE101089	AI010721	AA799666	BE108840	M94040	U21662	U25808

TABLE 3-continued

TABLE 3-continued

	121	DEL 3-con	illiucu			171	DLL 5-con	imucu	
	7478 genes	used to derive	RTI signature	<u>s</u>		7478 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #
AA944403	H35178	NM_012669	AA893505	AA924151	AW527421	AI412086	AW918940	AI105243	AA799709
AI704799	M11942	Y12517	AI058276	AB003400	NM_012776		BE110557	AI233763	AI070397
BF398009	M73808	AA819729	AI172267	AI227672	AI177863	BF285334	NM_012875		AI102943
D85435	Y12708	AF054826	AI177016	AI406500	AI406964	BF288060	AF095741	AA866432	AI231777
AI233765	AB017638	AI180337	AI233728	AW253963	AI411212	BF290997	AI231196	AA946017	BF551377
AA800539	AI169242	AI234533	AI406932	AW914642	AI556246	BF407158	AJ245646	AI105117	AB018791
AA892044	AI233232	BE105565	AI412180	AW918527	U62940	BF420447	AW525089	AI598410	AI008971
AA942808	AI237681	BF564263	AW143212	BE101505	AA800570	BF556463	AW528792	AW141364	BE102266
AA946508	AW143114	NM_012866		BF282984	AA946434	AA998047	BF410042	AW532663	BF399504
AW915559	AW913929	NM_019152		C06665	AI407954	AI231781	NM_017169		AA800001
BF282349	BE113268	AA944438	AI060197	U44979	AI170671	AI236726	AI227832	AA998468	AB010954
BF398332	NM_017048		AI230388	AA942726	AI409070	AJ003004	AI104378	AW142276	AF179370
AW143568	AA957492	AF110025	AI408502	AA944828	BE113315	AW531275	AI170657	AW914992	AI171990
AW916347	AF094609	AI145630	BE108850	AI169053	AA818128	AW918257	AI230061	BF285344	AW915681
BF408957	AI009654	AI176996	BE329450	AI171242	AB028626	BE108494	AW921738	BF561196	AW918311
AW915669	AI013906	AW141869	BF398626	AW915015 BE110412	AW915815	BE111850	BF419366	AB017188	NM_017182
AA800699 AI011749	AI171617 AW531909	M54926 AI169607	AA819234 AI103467	BE113269	BE101212 BF559919	BE113375 BE120015	AA892780 AA875425	AI406280 AW915764	X93352 AA924980
AI104431	BF393126	AI169746	AI103407 AI177412	U78977	AA875045	NM_017264		AA945568	AF172640
AI170825	X13058	AW915955	AI229902	AA848503	AI137420	AA926279	BF398144	AI176477	AI101380
AI575445	Y17326	BF282899	AW915152	AF244895	AW251313	AA946382	BE101784	AI599407	AI179992
AW251630	AA801094	BF400575	AW916942	AW435017	AW915638	AB008571	AI111559	BE113340	AI717425
BF287135	AI169140	U64030	AW917815	BE108968	AW917594	AI013657	AI169149	BF549893	AW916433
BF420680	AI232722	AF259504	AW919586	BF405135	BF524281	AI176468	AI175019	AA892273	BE098799
BF548086	AI236270	AI171230	BF291214	AA850872	BF556698	AW520324	AI177410	AA899959	BF397603
AA942949	AW915791	AI229647	AA943831	AA944332	AI233262	BF283406	BE107245	AF285103	AI102027
BF388434	U41803	AI235502	AF034582	AI600085	AI233718	BF418890	BE118650	AI176465	AI104258
AA892829	AA893241	AW523709	AF077000	AI600108	AI598371	BF420754	NM_012985	AI411365	AI454943
AB002151	AI228540	BE108860	AI412298	AW433865	AW141873	BE106191	BF407964	BE100986	AI059108
AI170414	AI317827	BF419854	L12384	AW913942	BE101766	BE099950	U61696	AA946375	BE110530
AI233729	AI575026	L25331	AI102688	AW916661	BF282301	BF407170	AA800277	AA955172	BF410389
AI236101	BE104107	AA818113	AI232248	AW921139	BF415017	AF110195	AA819086	AF255305	AI412276
AI412255	BF282890	AF056034	BE103304	BE101171	BF420639	AI012785	AI172459	AI169359	AW433846
BE101485 BE110671	BF287032 BF398047	AI407095 AW915655	BE109671 BE112899	BE106523 BE107223	NM_012789 NM_017299		BF397894 AI407555	AI408455 BF396218	AF002251 AI104146
BF283122	BF419646	BF387477	NM_017361		AI411436	AW 233983 AW 914085	AI556546	BF548597	AI104140 AI454536
BF414192	NM_017348		AA849734	BF563261	AJ303456	BE112582	AI577393	BF557304	AJ005424
NM_017013		BF555532	AA924717	AI406651	AF044058	BF286237	BF281749	U93692	AI233276
AI102745	AI105205	L37085	AI232347	BE096104	AI410001	BF399633	AF144701	AI232657	AI716471
NM_021676	AW918593	AA817752	AJ004912	BE101124	AW525660	AA945062	AW141326	BF414143	AI230758
AA799550	BE100453	AA858600	BE095620	BE109118	BE349725	AI103988	BF557792	NM_012960	AI412949
AB008538	BE102815	AI169490	BF398121	BF284014	BF282686	BE109599	BF420654	AA891821	AI600036
AF334379	BE103430	AI575402	BF417396	BF396629	BF549603	BF523605	AI059234	AW917596	AW253367
AI235934	BF282594	AW143173	NM_017015		BF407149	AI175803	AI232643	AI100850	BE104143
BF411031	AW919578	BE108396	AI103129	BF289928	AA924654	AI556502	BE113423	AI102689	AA799783
NM_013222		D21800	AI234816	BF565365	AW144382 AW915749	AI599995	AA892993	AI179136	AI716491
NM_019259 AI412015	BF283091	AA799499 AA892127	AI175507 BE119615	AI111991 BF286941	BF281388	AW917738 BF284345	X13817 AW915662	AW253642 BE118414	AW921162
AI169353	NM_019206		BF408841	AF200359	BF282084	M62388	AB006450	BF404027	AF110026 AI013011
AW252811	AA963094	AF311055	AI137756	AI009363	BF283385	AA924152	AI233857	BF414266	AI411227
NM_012619		AI169365	AW434991	AW915716	BF400719	AI600216	AW915056	AI412024	AI101580
NM_012019		AI407130	NM_019238		AI177621	AW523737	AI171211	AW919474	AI598381
AA851239	AW918097	AW527971	AF069306	BF523646	AI575104	NM_019144		AA801308	AW920761
AA899150	AW919037	AW916168	AI599945	AA894030	BE112007	Y00350	AF032120	AA818914	BF558116
AI171607	AW919937	NM_021745	AI137114	AI713140	AA848795	AA893208	AI169648	AF120111	AI555567
AI172029	BE118683	BE115558	AI232357	AW915146	AA894262	AI703715	AW918604	AI102947	BE099224
AI180458	AI598320	AA964789	AI412958	BF412293	AI230432	AW916925	BF397588	AI409731	BE112202
BE102485	BF281741	AI169729	AW251310	AB037424	AI548620	BE099060	NM_019213		BE117946
BF550566	BF285339	AI172272	BF417793	BE110618	AW917543	J05405	AA894297	AW913868	BF282388
BF556846	BF549027	AI179472	BF419240	NM_017326	BE115626	AA799331	BE104415	BF398537	AA800521
NM_013033		BF284775	U19614	AI073176	AI009222	AA944053	BF282678	AW526283	AA849788
AF030358	AA818203	BF398680	AW525945	AI411198	BE108018	AF184893	NM_019334		AF281304
AI176121	AW916939	BF410951	AA801212	BF398587	AI235192	AI172269	AI169328	AI412537	AI010455
AI598881	BE113338 BF408856	AI175767 AI599956	AI639285 AA800191	AA955157	BF283084 NM_012595	BE112892 BE410731	AI172092 AW528057	AI716902 AW434064	AI144663 AW915194
AW143543 AW915481	BF548630	A1599956 BE100802	AA800191 AA800535	AI105145 AI231011	AI178818	D50696	AF026476	AW 434064 BF396493	BF544320
BF399447	BF557395	BF407563	AW142925	AI231011 AI236640	AW525229	AB032178	AF136585	AA800576	AA944449
NM_017201	BF568009	AA893590	BE108810	AI412002	C06787	AI012381	AW918068	AI579376	AW142350
NM_017281	M29472	AA944576	BF399618	BE110561	D83948	AI180252		BE113316	AW531382
Z83044	U75928	AI169375	X67654	BE111986	Z71925	AI228249	AA925303	BF406661	AW915412
AI058938	AI009818	AW521376	AA893532	U15138	AA945915	AI230278	AI007987	AI233172	AY017337
AI137569	AI317880	AW918620	AA944158	BF397956	BF287826	AI408770	AI229046	AF110732	BE096311

TABLE 3-continued

TABLE 3-continued

	7478 genes	used to derive	RTI signature	<u>s</u>	7478 genes used to derive RTI signatures						
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #		
M75153	AI409748	BF420055	AI102991	BF417071	BE109616	BF410786	AF024622	AI169399	AF008554		
M83675	AW433870	AW143287	AA891940 BE100586	AF315374	BF419319	AI105417 V01224	BF405059 BF563404	BF281357 AW918387	D50864		
AA858879 AI231773	BF419628 M61142	AI105345 BF413977	AI233751	NM_017177 BF404344	AF136583 J03621	U07609	BF400779	BF408867	AF020346 NM_019145		
AI232273		BF398712	AW916097	AI555009	U02315	AI168935	NM_019314		NM_021766		
BE107540	U66322	AI408162	BE109950	AW919046	AW919325	BF281135	BF419671	U61184	BE101094		
BE113490	AI406508	AW523409	AI172301	BF549833	BF282951	L12025	NM_019179	AB006614	BE109569		
BE120629	AW915566	BF283600	AW520767	AA850288	NM_020080		BF397726	AW920729	BE117893		
L11004	BE115600	U69485	BE109512	AI411153	X98746	BF564840	NM_020301		L33413		
X74226	BE116507	BE109521	BF420279	AW916463	AF154914	AF163321	D38104	AA819339	AA818892		
AA858867 AA859922	AI171632 AI007841	AA944494 BF282132	BF393934	BF282695 AI410079	L27651	U38938 X06942	NM_019157 AI409500	AW918535 NM_012587	BF558524		
AF067728	AI599286	BF417400	AA800258 AI171764	AI411278	AA875041 AW919685	BF404901	AF100421	NM_012387			
AW920774	BE349648	NM_012891		M62763	D14437	D45920	NM_019290		AI175907		
AI412491	AW143711		NM_021264		X99338	NM_020074		AW251633	BE097102		
AW915621		AI406655	X52477	NM_013098	AF016180	AJ000555	AF141386	BF542467	BE111729		
BE101165	AI411113	U02096	NM_013104	BF396151	AI500969	AA859556	BF403998	BF565649	AI172498		
AI145899	AW913987	M31176	NM_020088		BE105541	U69550	NM_019272		AW915002		
AW917752	BE095840	U22830	AI009128	AF035963	BE108368	AI071605	U12402	BE120309	AW140991		
BE115557	BF411317	AW143269	NM_012629		U49235	BF557670	AI598429	BF388912	BE107195		
AA819400	BE101129	M55050	NM_013041		U66292	AF188699	BF414004	AW434670	BE117687		
AB049151	BE100823 BE101292	AI548036 M98820	J04731 NM_013178	AF054586	AI236780 AI599365	BF284311 AB000216	BF549324 Z49762	BE110658 NM_012900	U41164		
AI172464 AW141870	U53512	AI007936	BE111869	AB037937 NM_017280	NM_012896		Z49762 Z50051	AF062594	AI228240		
NM_019252		AW141286	AF012891	AI599294	AI176810	AW919217	U37026	BF420163	BF392884		
BE112384	AI410481	BF556841	BE109711	AI007877		NM_012526		U00964	BF546209		
BF285023	AB041723		NM_017170		BF415072	AW921292	AI235610	AI044845	AW918841		
AA800665	BF281200	AF176351	AI177747	AA859768	AB020019	NM_013078	NM_013224	BE113165	U25967		
AI178806	AA943011	AA943126	AI176502	AF168795	AI170357	BF555189	BE108347	BF393078	AI137259		
AI406906	BE096986	AW143102	AI105086	BF284341	AI716535	AI045026	AI137751	BF558742	AA944308		
AJ225623	AI044721	X78689	NM_012971		BF405610	U79031	AW917981	Z19087	AI407975		
AW918039	BE116383	BF550800	AW916756	M14952	AI176718	M60753	BE116152	AF150741	AF168362		
BF407819 AA849757	BE111699 AI104034	NM_013023 NM_012844		AI411426 BF285557	L27112 U08255	BE113295 AI111803	BE120545 AW919982	U57049 AA800382	AI406533 AW141129		
AI170714	AI548730	NM_013191		U34843	Y11490	AI111603 AI233752	BE109277	AW917673	AW915546		
BE109614	BE113022	M22926	BE111296	AF007789	BF408271	L07736	U93851	J04811	BF288240		
BE116918	BE113201	AI408780	M63991	D10693	D16829	BF401764	M87053	AW918559	U22424		
AI011510		NM_012822		D88672	AA946492	AF072835	BF283410	BE099796	AA800389		
BE115034	AI170354	U53449	U76551	AW918103	AI598346	X63995	M91597	AB046606	AJ223599		
NM_012670	BF411424	BF285022	D00569	AF003598	U28504	BF407531	AF030378	BE098930	NM_019147		
X96663	J03624	M96548	BF399655	NM_019241		NM_012699		AF079864	AI045083		
AI145851		NM_012521		NM_012694	AA946350	AA891949	M22253	BF557269	AJ132230		
BF547710	AW143082 AF199322	AI715955 BF404304	AI170387 NM_013154	AI233253	BE108246	NM_019371 BF417565		AA851305 AF201901	NM_017020		
AW528625 BE349838	Z46957	U21954	AW915339	NM_012702 NM_012716	AW141135 U81037	NM_012892	BE108756	AI411955	BE120578 NM_017010		
BF389726	L06040	Z96106	AW919159	NM_019223	AA946467	BE107187	BE109242	AW143142	BF282689		
BF523622	AF180350	BF550451	AI172174	Y07704	AI412189	BE108224	AW533482	BE109664	D17309		
BE111787	AB040802	AF135115	U04998	X06423	AI180349	L01702	BE111827	AA900180	L36088		
AI170768	BE109138	X92097	AF069770	Z18877	X95189	AI177168	NM_017149	AF230638	NM_012700		
BE113043	NM_021695		D12978	AI170265	BE117941	M59967	AI598306	BE101140	NM_021593		
BF282314	NM_019125		BF288153	U28356	AW916860	X98517	BE103482	D37979	U48245		
M57547		AW915423	AI172352	AA945099	BF523660	AJ002556	BE104535		NM_017154		
AA858518	AW141928 X00469	AA893251 AI229720	NM_012744 U57063	X95096 J04628		NM_013172		AI711110 NM 017286	U89695		
AI575433 AA818520	AA799329	D49494	NM_019189		U92802 BF285568	AA943114	X13722	NM_017286 NM_019137			
AA893517	AA818947	U44125	L26009	NM_013100	BF283308 BF281914	BF401491	AI178476	AI411375	J03025		
AF165892	AI178768	X67859	BE096501	D63834	U59672	BF287843	U56241	AI236120	AA946014		
AI179365	AI010317	BF394161	AW918276	AW920575	U90829	BE095997	U67137	AW915825	BF287827		
AI230346		NM_021669		AF203374	AA858794	BE109638	BE108748	BE113365	BE098800		
BF406604	BF563201	AI406856	NM_012826		AI317854	L14851	AI172248	BF393902	AI712686		
BE113454	D10041	BE107032	BF420059	M31155	AJ009698	M83679	AW526136	BF556880	BF404426		
AI171781	AF097723	NM_017158		BE113362	NM_012742		BE107192	AA894099	J02962		
AI179316	U89744		NM_017076		AB019120	BF399595	U17253	BF417476	NM_020471		
AI171367	M58364	BF549748	AF013598	AW920993 AI176592	AW918468 J04147	BE112983	X95577	U53882	X89383 U25651		
BE109901 BE329347	BF398051 AW434139	NM_017136 AW143169	AF242391 AI170665	X54467	NM_017129	AI232716 BE100617	AI600081 AA819316	AF153012 BF288244	AB015308		
AI410096	NM_021656		AB018049	AF150106	U25055	BE100017	BE113655	NM_012773			
AI411531	AF205717	Z14119	AA801173	AJ002745	NM_017278		AI233213	M30596	BF284768		
BE110545	L33916	BE113272	AA818949	BE104375	BF398182	BF399649	AI178556	AW142966	BF406261		
BE111677	AW527564	BE121346	AJ132352	D83792	AW915563	NM_012707		AI228955	D30666		
AW141664	L27059	D14048	AA996961	AF082534	AA851914	BF398696	BF408425	AW523755	Y00697		
AF000973	BF289566	U31203	AA800501	BE118055	U56859	M88469	AW142667	AW918188	AB022714		

TABLE 3-continued

TABLE 3-continued

	7478 genes	used to derive	RTI signature	<u>s</u>	7478 genes used to derive RTI signatures						
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #		
AA848534	AI236753	BF404316	BF522212	AF132046	AI406532	U30831	NM_012807		BE113015		
AA944398 AF022247	M83196 AF067793	BF555947 NM_019165	AF036335	AI178938 M22923	NM_017110 NM_021774		BE104961 BF281284	BE109712 AF254800	BF283407		
AF022247 AW434092	AI137506	U92072	AF083269	NM_012591	AI103955	AI407187	AW915886	AI406670	AI639411 AW526270		
BE108809	U48592	AF156878	AI235493	AF093536	AA818197	BF408452	BE120498	U20999	U08257		
BF404853	BE113616	BE114418	AI411056	AI406525	AW520354	BF413245	AI175533	AB006137	AI716159		
AB021971	X07467	AI179460	L27058	AI408017	Z21513	M88096	NM_020308		AW143164		
BE113076	BF285915	NM_020081	AW253040	BF284776	BE097309	NM_017179	AI145869	BF400588	AW919130		
BF414136	BF563467	AI102524	AA818020	BE113205	BE118454	AF279918	AW142932	BE097982	BE107279		
AI145380	D16465	NM_017150		M14050	BF290106	AI007974	AW143294	M85183	AI236615		
AA943794	NM_019318 U90888	AB022209 BF285150	AW527440 BE113660	AA892918 NM_013050	M31837	AW144302 BF402472	AW251657 AW525099	AA892987	BF291260		
AI146056 AI178808	AA899951	AI411991	BF389120	U42413	AA955605 AW918716	BF416877	AW523504	AF106325 BE117939	NM_019350 AA858745		
BE109381	BF556350	NM_012734		AB016160	BF288254	AI411670	AW915685	BE119991	AI598946		
M55534	U17971	AW142654	AW918441	AI410127	NM_012811		BE115417	BE126380	AJ245648		
BE329415		NM_012610		NM_019384	AB010960	NM_017323		AB033418	AW530332		
AB003478	AA874975	AI176548	BE109561	AI102061	AF081582	X73653	AI177360	AF068861	AA957010		
AF157511	AF110023	AA850242	NM_017003	AA946349	BF393949	AW251852	NM_021586		AW915165		
BF408990	AJ225654	BF396462	AI236928	NM_019256	AI009820	BF281178	AI598507	AW918920	BF564460		
BF551318	AW915004	NM_012913		AF008114	AI229209	AJ002942	AW915843	X73292	AA849731		
L05084	BE097085	AI145761	AI407992	AI230591		NM_017214		AW143907	BF543359		
NM_021741 BF553500	AA801331	AI411297 NM 017060	Y09164	BE101290 NM_021661	NM_016994 AW533098	NM_017035	AI176646	BE113277 BF389143	AA944327 AW918368		
BF564759	BE109744	AF281635	X05341	NM_017062	BF404778	AF100172	AI409051 AI409861	AI105450	BE100015		
AI176478	AI713217	U42388	AI411422	AI408969	AF032925	AW142717	BF281872	AW526673	BE100965		
AI454928	BF396314	BF409208	AW919170	AW918198	AA800046	AI409727	AF026505	AA848804	BE110621		
AI599484	NM_012921		AI763826	AA800744	NM_012653		BE103793	AI411527	BE112971		
AW917650	U89280	AW143820	BF288140	AB028461	NM_013151		AI169749	AI235294	AI102009		
AB001089	AW917574	BF282574	X96589	AW143077	AI176727	AW140397	AI454845	U66461	AW915541		
BE098025	BE105864	AF178689	BF406991	BE101126	AW143676	BE108985	M25073	AA945069	AW921544		
BE101151	U14914	BE107098	AA817836	BF555858	BE109246	BF285313	AW916692	AF081503	BE100576		
BF282700	AI227916	BF407134	X15834	BF556693	AF176072	BF289100	AI177083	AI408928	BE113248		
M74067	AW920324	L39018	BF555544	AI717140	AI170289	D10554	AI228159	AW920454	Z29486		
AI232138 BE110691	BF281577 AI600068	BE097840 BE107410	AB046544 L13041	AF240784 AW916911	U27186 AW530292	M83107 AI411222	AW434103 BE108388	BF414412 BF419792	AI169653 AI172320		
BF282674	BF282471	AI227686	AF009511	BF284509	BF387258	BE107075	BF282088	AW252109	AI76972		
AF015949	BF396350	AW916943	BF567426	BF418775	AI175556	AF218826	BF420183	NM_012904			
AF054870	M91214	AF227741	NM_017161		AW143156	AI102758	AI409150	NM_021757			
AI009608	X82021	NM_021670		AA800737	AW917977	AI231782	AW915035	AW916774	AW918991		
AI411793	AF062389	X52196	AF063851	AI146063	X14773	D88190	BF284983	AW915540	BF551808		
	BE116153	AB025784	BF415013	AI407061	AI230732	AI555819	AW435159	BE096098	U36482		
BF282483	AA851302	NM_021776		BE115604	AI412740	D21158	Y00047	AI104256	AA997745		
AA800364	AW914045	U03491	AI233343	BE116180	BF406286	M77246	AI411360	AI179391	AI069922		
AA945579 BF392911	BE113101 BF407194	X02610 AA818377	BE109221 D78610	BF419010 BF557276	NM_020101 U81160	AII /1651 AI548655	X78855 AW251401	AI556315 BE109108	BE107622 BF567869		
NM_019122		AI599232	AW917185	NM_012761	AI113186	AI555457	Z49858	BF417010	AW919868		
AA818602	AI703713	AW915797	NM_019368		AA893584	BF404590	AA944079	BF548406	BE118251		
AI172262	BE109919	BF285528	AA800815	AA818342	AF161588	BF411162	AI172266		AW915682		
M85299	BE096027	BF566689	AF061873	AI230596		NM_016996	BE107489	BF417442	AI136709		
NM_012770	AI409037	AI574745	AW918233	AI406712	AI170410	AI101373	BE116220	AI407858	BE110614		
AA955527	AW143190	BE112950	BF542548	AI410452	BE096257	AW918990	AA945713	BF282381	AA899160		
AF222712	BF396115	X65948	AA851282	AI548615	M29294	BE107805	AW920687	D29960	AB029559		
AA943573	Y15054	AW252115	AW918031	AW251686	U36992	AI177022	BF415222	L26267	AI409045		
AI412866 BF388797	AA946222 AI511282	BE097244 NM_021590	BE096652	NM_017300 M64381	X66022 AA848437	AI411391 AA891551	AI177058 AW532870	AI229821 AI230362	AF117330 BF523591		
BF403842	AW915567	BF282395	NM_021866		AA850317	BE106307	BE108835	AA891733	AI411845		
NM_012624		BE116554	AW527690	AW919837	D13518	BF393917	NM_020075		AI502504		
U78090	BE115948	U55836	AI180050	NM 017238	AI008964	BE110638	AI102743	BF281852	AW143215		
AI176298	AW915148	X16262	AW917831	AI012336	AI639504	M69246	AI169386	AA 799700	BE098555		
AW523419	BE107295	BF550679	BE109095	AW919892	AW919920	AI227943	BE109681	BF414012	AI009644		
AW917114	BF395101	AI170390	BF408792	NM_017243	AW920557	BF524978	AI599479	BF393863	D87950		
M22756	AI146156	AW915776	AF221622	AF194371	AI008371	U31352	BF283237	AB018253	AA893640		
U48702	AI599410	BF284067	BF388220	AI102802	AI230578	AI010423	NM_012797		AI009197		
	AW144095	NM_019376		AI138048	AW915692	M13979	AI171794	M34477	AI171088		
AA944552	AW915236	AI072251	AI178019	L38644	BE099875	AI177590	U27191	AI236691	AW921523		
AI408249 U25281	NM_012969 AW253398	AI007768	BE111361 BE111820	NM_013124 X54862	BF282648 AA850576	AI111863 AI235282	L27339 AW916023	AI406487 BF283073	BE113399 AF016252		
NM 012889	AW 233398 BF282987	AW921975	AW915041	AJ4802 AI178452	AW915782	AW919696	BE098778	AI411772	BF420684		
AB009372	AF234260	AA925469	AW915273	AI578745	BE102100	BF553948	AF136943	AW251612	BF557013		
J00750	BF416935	AI102804	AI556256	L22294	D17711	AA858925	AI178272	AW525370	AI103146		
NM_012834		AW141615	AA817802	AI178361	AA998252	AF033027	AI231505	BF551370	AI501407		
AA848305	BF553981	BF404819	AI176838	AW520781	AW531386	J05035	BE096516	BF393595	AW528778		

TABLE 3-continued

TABLE 3-continued

	7478 genes	used to derive	RTI signature	<u>s</u>		7478 genes	used to derive	RTI signature	<u>s</u>
Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #	Accession #
BF402375	AI710879	NM_012853		NM_017282	BF411461	BF397012	D14046	AJ238278	AW916153
AI101189	NM_019163		L28801 NM_017037	AA893192	AA963071	BF550545	X79807 AW530379	J04487 AF008197	AJ301677
AI410802 AI599568	U71294 AA849987	AA944485 U26595	BF285026	BE107103	AI171951 AI410391	D10854 AA818089		BF396279	BF566748 NM_017139
L14684	BE101435	AA818910	AA849958	BF415001	NM_012552		X68282	M59814	BF406213
AF324255	NM_020076		BF550748	NM_017209	AI009200	Y07744	AF062402	AI230548	D10233
U82591	BF283735	AB048711	AI409857	BE111727	BF281133	BE114123	AW251683	L20823	BE111690
X59601	BF394563	BE121333	BF394528	BE116914	NM_021752		NM_012852		BF410771
AF214568	AW917562	AJ225647	AI235367	BF398071	AW915445	BF289154	U12571	AW143086	AI045074
AI406304	AW918237	BF283417	BE104454	AI007924	BF284328	BE111731	NM_017304		AI137283
U57715 BE097153	NM_012827 U73503	AF216807	AW915120 BE117683	AI233773 BF281865	BF285068 BF288092	BE119400 BF393862	NM_019175 AI013038	U23438	NM_017049 AA891922
BF281969	AI170258	AF323174	AF025424	AI175544	AW918538	AW251199	X70706	NM_021701	
AI716436	AW914009	AA849774	AA859141	AW914973	BF282327	AW526089	AB032827	AW534166	AI170382
BF417252	AW914041	AF003926	BE116512	AW916823	BE101579	AW919172	BF286192	L15619	BF398540
AI010312	BE109628	BF557691	AA946032	BE099999	BE116560	BE110609	AA892531	BF398602	AI145586
BE111694	BF555169	NM_013027		BE110645	NM_012758		U52103	X15800	BE110674
BE113210	AI179795	AA848530	BE108886	AF002705	AA996543	AW434308	BF412297	M81642	BF419044
BE117891	BE102427	AW529231	BF551148	AW918614	BF415031	BF554891	AI144958	NM_012790	
AW435036 AA943752	BE111811 NM_017313	BF557674	AA942690 AW914062	AW917503 BE113989	AI410349 AF134054	AW 254375 AI 235510	AF291437 U73174	AB030947 L08814	BF283759 L11007
AA800739	AA819318	AW919277	BF408129	BF284855	AW916920	AF219904	AI176327	X58828	NM_017248
AI070523	AW524523	AW918452	BF550572	AW915142	AI385370	BE115570	U06273	NM_017017	
AI177379	BF411622	BE100774	AW253429	AW915484	BE100607	AI234173	AF095449	AF035156	AB001321
AI406469	BF558866	BF396644	AF223951	AI411205	AI233786	AW525042	AF269251	AI715321	AI406494
AI598307	U57391	AA892362	AF277899	AF036255	AI713324	BE104891	NM_013127		BE097282
L19699	X62277	AI227985	BF557299	AI579643	BF290678	X71429	X97831	X52590	AW435315
AA901066	AA892240	BE108201	AW915121	AA799544	BF414193	BF283351	X53477	BE102840 NM 017333	D37920
BF412565 M84719	AW251483 AW525372	AA799313 AI236027	NM_017349 AI406938	AI102643	AI406350 AI411530	BF415054 BF282715	NM_017284 NM_013225	_	BF400782
BE110574	X67156	BE120602	AI010033	AF302085	BE111849	AI010351	U66566	AF011788	AI231089
BE120360	AI231309	X04070	AI407930	AW251641	AW916376	AI112973	S79760	AB032551	BF417360
BF282675	BF420064	X54640	AI599023	AA943815	AA799532	BF558479	AW143231	BF549260	NM_017173
AI408993	D50695	AA894305	NM_012680	AW253646	AI170763	X53232	BF564152	U53706	X98399
NM_019201	U96638	AI229183	Y00826	AF221952	AW141730	AW144324	AW915661	X68191	Y13588
AW144669		AI411790	AW254166	BE109709	BE111512			M14053	BF283107
BE110128 BE113321	AF274057 AI409506	AW143855 BE118562	BF411147 AI408520	BF549441 AI172472	AI102788 BF284711	AI169374 AI555565	BE115943 AW920769	NM_019352 BF403999	AA/99338 BE116101
BF392443	AI409300 AI412429	AA892772	BF283250	NM_013113	AI103993	AA800637	BE110731	NM_012930	
BF556332	BE116848	AA892922	AW917768	AI228624	BF557930	AI009796	Y16641	AW920609	AF016387
X66370	BF551342	AA800249	BF396485	BE112237	BF563406	AI407449	M87067	AW918854	AW918052
AI178214	AI137218	AA858572	BF401591	BF281954	BF288328	BF287028	AF016183	NM_013174	L35771
AI230884	AI171769	AW918182	AI169706	AF170253	L26288	AI145385	D86373	AA964289	AF007108
AI556402	AI412763	D84667	BE109678	AI012264	M34043	H35156	J05181	AB015746	AW142311
U35245	AA818692	U89282	NM_017163		AB046442	BF284885	BE110547	U78517	AA819501
U40188 AI169635	AI178158 AI231286	AI010433 BE113013	AA799691 BE104865	D38082 AA946389	AF239045 AI639012	AI010413	NM_021696 NM_021758		AW144294 Y00752
AI410456	AB006914	BF396371	BF284093	BF523723	AA942681	AI175064	NM_017230		NM_012651
BF288651	AI411156	D13124	AI013075	M57728	AF276774	BE111765	NM_017127		AW251942
BF397951	AW916461	AA996628	AW433875	AA818582	BE099976	V01222	BF557572	U93306	BF399135
AF000578	BF403712	AB047556	AW433883	BE097279	AI412079	BE116947	U05989	AJ010750	BF289492
AW520760	AI102486	AI137161	BF418588	BF550623		NM_012883		BF281419	AA946518
BE104290	AI137233	AW252891	AA955616	BF407480	M75148	NM_017094		AI715893	BE107610
BE117164 BF406590	AI175494 BE109633	BE111879 BE115051	AF227200 BE097298	AF087431 AI232979	NM_012594 AW143162	AW143513	AI406390 BF407501	BF407203 AI237636	
AI229833	AB000199	AA945898	AW914090	AW251107	BF566346	AF187814	BE109573	AF095740	
AI009089	AI009094	AI012613	M61219	AW253004	BE101876	U51583	AI009156	AI179711	
AI012598		AI407067	BF406413	BF419187	BF551593	BF567845	AW917598	AW527815	
AI228598	AI408482	BE116889	AW143981	AA891860	AF087433	AW920343	BF289001	AA945149	
BF399993	Z11994	X71873	AW919527	AI172500	AA943764	BE105589	BF281975	AF234765	
L12382	AI178922	AA900562	AA893621	AA945753	NM_012747		D83538	BE110624	
AI176042	AI236063	AI170409	AW141940	AW252110	AI179315	AW913871	AI177053	BF406562	
BF562819	AW526015	BF400995	AI137912	AW254010	BF548006	AI009007	BF393285	D00859	
X52140 AA996888	BE109952 D26180	M27893 BF284313	BF406514 U00926	AW434299 M23674	AI175454 AF013967	D21132 AW253843	D82928 M29293	BE109704 AI407113	
AI105088	AA851369	AI230697	AI178155	AI169176	L22022	BE109532	BF291213	BF396256	
AI172150	AI227996	AI575056	AI716607	AA851241	NM_017190		AF017756	AI180187	
BE109232	AA858649	AW914867	BG153368	BE116927	AI575072	AW919439	BE118425	BE109634	
BF420717	BE096995	M33648	AF269283	AI548722	BE111659	BE098855	BF556755	AA944176	
BF282933	AI716456	AB032243	NM_021678		BF283830	AA848420	BF282147	BF395125	
AI228642	AA957770	L07578	X68101	BE109116	BE109642	AI406499	BE108922		
AI599819	AA945696	NM_017210	BF567821	AW522132	AI070732	AI406520	BF402664		

TABLE 3-continued

Accession # Accession # Accession # Accession # BF283754 BF410170 L22339 BF405725 AI232332 NM_013177 NM_013186 U67140 AF110024 BF420043 AA893191 AW143526 M55250 AA817813 BF555225 AA799784 BE111345 X71916 A1013110 AA817817 A1070303 BF412643 Al179413 AA965185 A1170664 A1231827 BE109656 M94043 A1579555 NM_017026 A1406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BE128566 AW141446 NM_019234 BF553144 AW915616 A1598719 AI231290 BE116574 AA801133 BF414497 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 A140827 BF548743 </th <th colspan="11">7478 genes used to derive RTI signatures</th>	7478 genes used to derive RTI signatures										
BF405725 AI232332 NM_013177 NM_013186 U67140 AF110024 BF420043 AA893191 AW143526 M55250 AA817813 BF555225 AA799784 BE111345 X71916 AI013110 AA817817 AI070303 BF412643 AI179413 AA965185 AI170664 AI231827 BE109656 M94043 AI579555 NM_017026 AI406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA8848639 M35495 L39991 BF5553984 BF411842 AB028933 <td< td=""><td>Accession #</td><td>Accession #</td><td>Accession #</td><td>Accession #</td><td>Accession #</td></td<>	Accession #	Accession #	Accession #	Accession #	Accession #						
NM_013186 U67140 AF110024 BF420043 AA893191 AW143526 M55250 AA817813 BF555225 AA799784 BE111345 X71916 AI013110 AA817817 AI070303 BF412643 AI1 79413 AA965185 AI170664 AI231827 BE109656 M94043 AI579555 NM_017026 AI406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AIS98719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA88639 M35495 L39991 BF395067 L14936 AA81768 AW918011 AI408827 BF548743 AI072	BF283754	BF410170	L22339								
BF420043 AA893191 AW143526 M55250 AA817813 BF555225 AA799784 BE111345 X71916 AI013110 AA817817 AI070303 BF412643 AI179413 AA965185 AI170664 AI231827 BE109656 M94043 A1579555 NM_017026 AI406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA84639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 A1716480 BF55	BF405725	AI232332	NM_013177								
M55250 AA817813 BF555225 AA799784 BE111345 X71916 AI013110 AA817817 A1070303 BF412643 AI179413 AA965185 AI170664 AI231827 BE109656 M94043 AI579555 NM_017026 AI406275 U46149 D89375 AF039203 BF309070 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 AI716480 BF5553984 BF411842 AB028933 AI	NM_013186	U67140	AF110024								
AA799784 BE111345 X71916 AI013110 AA817817 A1070303 BF412643 AI179413 AA965185 AI170664 AI231827 BE109656 M94043 AI579555 NM_017026 AI406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 A1045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 AI716480 BF5553984 BF411842 AB028933 AI235222 AI178134 AA859631	BF420043	AA893191	AW143526								
AI013110 AA817817 AI070303 BF412643 AI179413 AA965185 AI170664 AI231827 BE109656 M94043 AI579555 NM_017026 AI406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA88639 M35495 L39991 BF395067 L14936 AA81768 AW918011 AI08827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 A1235222 AI178134 AA859631 VM_017106 BE107324 D85189 U78	M55250	AA817813	BF555225								
BF412643 AI179413 AA965185 AI170664 AI231827 BE109656 M94043 AI579555 NM_017026 AI406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF444997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI2355222 AI178134 AA859631 VM_017104 BE111770 AI411088 AA900434 AA8800708 AI407320 AF049344	AA799784	BE111345	X71916								
AI170664 AI231827 BE109656 M94043 AI579555 NM_017026 AI406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BEI11770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 A1007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AI013110	AA817817	AI070303								
M94043 Al579555 NM_017026 Al406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 Al045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AL598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 Al170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 Al408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI233452 AI170376 A	BF412643	AI179413	AA965185								
AI406275 U46149 D89375 AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA88639 M55495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 VM_017160 BE107324 D85189 U78875 BEI11609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI1	AI170664	AI231827	BE109656								
AF039203 BF390970 BE100771 AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 A1408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 VM_017104 BE10170 BE107324 D85189 U78875 BE111609 NM_017104 BE11170 AI411088 AA900434 AA880708 A1407320 AF049344 M27905 A1233452 AI170376 AI172075 AA880487 A1007704 <	M94043	AI579555	NM_017026								
AA818364 BF405581 U54807 BF393486 AI045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AI598719 AI231290 BEI16574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 A1716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BEI11770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA	AI406275	U46149	D89375								
BF393486 Al045035 X99326 BE128566 AW141446 NM_019234 BF563114 AW915616 AL598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 AI716480 BF5553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 A1007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 <td< td=""><td>AF039203</td><td>BF390970</td><td>BE100771</td><td></td><td></td></td<>	AF039203	BF390970	BE100771								
BE128566 AW141446 NM_019234 BF563114 AW915616 AL598719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE106609 Al170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 A1408827 BF548743 AI072218 AI410818 AI716480 BF5553984 BF411842 AB028933 AI235222 A117814 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA893164 AS1070 D50694 AF051943	AA818364	BF405581	U54807								
BF563114 AW915616 AI558719 AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 VM_017106 BE107324 D85189 U78875 BE111609 NM_017104 BE11170 AI411088 AA900434 AA880708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA893164 AA8000707 D50694 AF051943 A	BF393486	AI045035	X99326								
AI231290 BE116574 AA801133 BF414997 NM_019124 U10894 AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 A716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BEI11770 A1411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	BE128566	AW141446									
BF414997 NM_019124 U10894 AW921546 BE100609 A1170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 A716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D88189 U78875 BE111609 NM_017104 BEI11770 A1411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	BF563114	AW915616	AI598719								
AW921546 BE100609 AI170303 AI044124 NM_013130 NM_019281 AA848639 M35495 L3991 BF395067 L14936 AA817968 AW918011 A1408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AI231290										
AI044124 NM_013130 NM_019281 AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 A108827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 A1178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	BF414997	NM_019124	U10894								
AA848639 M35495 L39991 BF395067 L14936 AA817968 AW918011 A1408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 A1411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800004 BF397919 AA8943164 AA800004 BF397919 AA894316 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AW921546	BE100609	AI170303								
BF395067 L14936 AA817968 AW918011 AI408827 BF548743 AI072218 AI410818 A1716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BEI11770 A1411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AI044124	NM_013130	NM_019281								
AW918011 AI408827 BF548743 AI072218 AI410818 AI716480 BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 A1007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AA848639	M35495	L39991								
AI072218 AI410818 AI716480 BF555984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	BF395067	L14936									
BF553984 BF411842 AB028933 AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AW918011	AI408827	BF548743								
AI235222 AI178134 AA859631 NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AI072218	AI410818	AI716480								
NM_017160 BE107324 D85189 U78875 BE111609 NM_017104 BE111770 Al411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	BF553984	BF411842	AB028933								
U78875 BE111609 NM_017104 BE111770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780											
BE111770 AI411088 AA900434 AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	NM_017160	BE107324	D85189								
AA800708 AI407320 AF049344 M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	U78875	BE111609	NM_017104								
M27905 AI233452 AI170376 AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	BE111770	AI411088									
AI172075 AA850487 AJ007704 M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AA800708	AI407320	AF049344								
M55075 BF283861 Y13380 AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	M27905	AI233452	AI170376								
AA944549 AA800291 AA893164 AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780	AI172075	AA850487	AJ007704								
AA800004 BF397919 AA894306 X51707 D50694 AF051943 AI179640 AI412931 BF558780		BF283861									
X51707 D50694 AF051943 AI179640 AI412931 BF558780	AA944549	AA800291	AA893164								
AI179640 AI412931 BF558780	AA800004	BF397919	AA894306								
	X51707	D50694	AF051943								
X74125 NM_017355 X61677											
	X74125	NM_017355	X61677								

[0142] The signature used to predict the presence or absence of future renal tubular injury was derived using a robust linear programming support vector machine (SVM) algorithm as previously described (see e.g., El Ghaoui, L., G. R. G. Lanckriet, and G. Natsoulis, 2003, "Robust classifiers with interval data" Report # UCB/CSD-03-1279. Computer Science Division (EECS), University of California, Berkeley, Calif.; and U.S. provisional applications US Ser. No. 60/495,975, filed Aug. 13, 2003 and U.S. Ser. No. 60/495,081, filed Aug. 13, 2003, each of which is hereby incorporated by reference herein). Briefly, the SVM algorithm finds an optimal linear combination of variables (i.e., gene expression measurements) that best separate the two classes of experiments in m dimensional space, where m is equal to 7479. The general form of this linear-discriminant based classifier is defined by n variables: $x_1, x_2, \dots x_n$ and n associated constants (i.e., weights): $a_1, a_2, \dots a_n$, such that:

$$S = \sum_{i=1}^{n} a_i x_i - b$$

where S is the scalar product and b is the bias term. Evaluation of S for a test experiment across the n genes in the signature determines what side of the hyperplane in m dimensional space the test experiment lies, and thus the result of the classification. Experiments with scalar products greater than 0 are considered positive for sub-chronic nephrotoxicity.

[0143] Signature Validation

[0144] Cross-validation provides a reasonable approximation of the estimated performance on independent test samples. The signature was trained and validated using a split sample cross validation procedure. Within each partition of the data set, 80% of the positives and 20% of the negatives were randomly selected and used as a training set to derive a unique signature, which was subsequently used to classify the remaining test cases of known label. This process was repeated 40 times, and the overall performance of the signature was measured as the percent true positive and true negative rate averaged over the 40 partitions of the data set, which is equivalent to testing 392 samples. Splitting the dataset by other fractions or by leave-one-out cross validation gave similar performance estimates.

[0145] Cross validation using 40 random iterative splits (80:20 training:test) resulted in an estimated sensitivity, or true positive rate, of 83.3%, and a specificity, or true negative rate, of 94.0%. Leave-one-out cross-validation produced similar results.

[0146] To test whether the algorithm is identifying a true pattern in the training set, but not a random data set, the labels for the 64 experiments were randomly assigned and a signature was derived and subject to cross-validation as above. This process was repeated 99 times. As expected, the average test log odds closely centered about zero (-0.004±0.86), with a range of -2.3 to 2.9. By comparison, the true label set had a log odds ratio of 4.4, which was significantly greater than expected by chance (p<0.0001).

[0147] Results

[0148] Using 7478 pre-selected genes whose accession numbers are listed in Table 3, the SVM algorithm was trained to produce a gene signature for renal tubule injury comprising 35 genes, their associated weights and a bias term that perfectly classified the training set. The 35 genes and the parameters of the signature are depicted in FIG. 1. Average impact represents the contribution of each gene towards the scalar product, and is calculated as the product of the average log₁₀ ratio and the weight calculated across the 15 nephrotoxicants in the positive class listed in Table 2.

[0149] As shown in FIG. 1, the genes are ranked in descending order of percent contribution, which is calculated as the fraction of the average positive impact each gene in the positive training class has relative to the sum of all positive impacts. Genes with a negative average impact are considered penalty genes. The expression \log_{10} ratio of each gene was plotted in the depicted "heat map" across all 15 treatments in the training set. The sum of the impact across all 35 genes for each treatment, and the resulting scalar product are presented along the two rows below the plot. The bias term for the 35 gene signature was 0.58.

[0150] The 35 genes identified represent 35 unique Unigene clusters. This 35 gene signature identifies compound treatments that are predicted to cause future renal tubular injury in the rat based on kidney expression data from short term (<=5 days) in vivo studies.

US 2006/0199205 A1 Sep. 7, 2006

[0151] The product of the weight and the average log₁₀ ratio across the 15 positive experiments in the training set indicated that 31 of the 35 genes are considered "reward" genes, as they represent expression changes that positively contribute to the signature score (i.e., the scalar product). The reward genes assure sensitivity of the signature by rewarding expression changes consistent with nephrotoxicity. A positive scalar product indicates the experiment is predicted to be positive for future renal tubular injury, while a negative scalar product indicates the experiment is negative for future renal tubular injury. The remaining 4 genes in the signature are considered "penalty" genes as they represent expression changes that negatively contribute to a scalar product. Penalty genes assure specificity of the signature by penalizing expression changes not consistent with nephrotoxicity.

[0152] The genes and bias term in the signature are weighted such that the classification threshold (i.e., zero) is equidistant, by one unit, between the positive class and negative class experiments in the training set.

[0153] Of the 31 reward genes, 15 have an average expression \log_{10} ratio greater than zero and are therefore induced on average by the nephrotoxicants, while the remaining 16 are on averaged repressed by the nephrotoxicants. Examination of the expression changes across the 15 nephrotoxicants in the training set reveals that most genes are not consistently altered in the same direction by all treatments (FIG. 1). Instead, it is the sum of the product of the weight and \log_{10} ratio (i.e., impact) across all 35 signature genes, less the bias, that results in an accurate classification. For example, Cyclin-dependent kinase inhibitor 1A (U24174) or the EST AW143082 are induced and repressed to varying degrees by compounds in the positive class, thus indicating that individual genes would be poor classifiers when used individually. This highlights the limitations of using single

genes for classification and also illustrates the basis for signature robustness since classification decisions are not dependent on any one gene that may be subject to experimental error.

Example 4

Stripping of Renal Tubule Injury Signatures to Produce a Necessary Set of Genes

[0154] In order to understand the biological basis of classification and provide a subset of genes useful in alternative signatures for renal tubule injury, an iterative approach was taken in order to identify all the genes that are necessary and sufficient to classify the training set.

[0155] Starting with the 7478 pre-selected genes on the Codelink RU1 microarray, a signature was generated with the SVM algorithm and cross-validated using multiple random partitions (80% training: 20% test) of the data set. The 35 genes identified previously in the first signature (i.e., "iteration 1" in Table 4) as being sufficient to classify the training set were removed and the algorithm repeated to identify additional genes. This identified an additional 37 genes (i.e., the genes in "iteration 2" in Table 4) that were able to classify the training set with a log odds of 3.80. This approach was repeated until the test LOR of the model reached zero, which occurred after 14 iterations and which consumed 622 genes. Based on the first 5 iterations, 186 genes were identified to be necessary to classify the training set with a test LOR of 1.64 (Table 4), which is approximately 2 standard deviations greater than the average LOR achieved with random label sets. Importantly though, it identifies a reasonable number of genes with a demonstrated ability to uniquely discriminate nephrotoxicants with an approximate accuracy of 76%. These genes are listed in Table 4.

TABLE 4

186 genes identified to be necessary and sufficient to classify the training set.										
Probe	Iteration	Weight	Impact	Mean Logratio Positive Class	Mean Logratio Negative Class	Unigene ID	UniGene Description			
AI105417	1	-0.89	0.261	-0.294	-0.172	Rn.8180	neuronal regeneration related protein			
BF404557	1	-1.36	0.213	-0.156	0.077	Rn.50972	ESTs			
U08257	1	0.88	0.149	0.170	0.029	Rn.10049	Glutamate receptor, ionotropic, kainate 4			
BF285022	1	1.46	0.143	0.097	-0.013	Rn.24387	ESTs			
AF155910	1	0.55	0.125	0.226	0.002	Rn.92316	heat shock 27 kD protein family, member 7 (cardiovascular)			
AI144646	1	0.63	0.108	0.171	-0.075	Rn.36522	gap junction protein, alpha 12, 47 kDa (Hs.) (DBSS_strong)			
AI105049	1	0.82	0.104	0.126	-0.018	Rn.23565	ESTs			
AI227912	1	0.46	0.074	0.160	-0.026	Rn.873	Sorting nexin 3 (SDP3 protein) (Hs.) (DBSS_strong)			
AW916023	1	-0.64	0.074	-0.116	-0.011	Rn.6788	Kelch-like ECH-associated protein 1 (Cytosolic inhibitor of Nrf2) (INrf2) (Rn.) (DBSS_weak)			
BF403410	1	0.42	0.068	0.163	0.020	Rn.23087	Homo sapiens clone 25048 mRNA sequence (Hs.) (DBSS)			
Y00697	1	0.63	0.067	0.106	0.048	Rn.1294	Cathepsin L			
AW143082	1	-0.30	0.056	-0.186	0.361	Rn.22057	ESTs			

TABLE 4-continued

Probe Iteration Weight Impact Class Class Unigene Unigene Class ID UniGene Description	
Alio2732 1 -0.31 0.035 -0.113 0.064 Rn.7539 ESTs Ali76933 1 0.46 0.035 0.076 -0.048 Rn.23658 ajuba (Mm.) (DBSS) AF208288 1 -0.27 0.034 -0.127 0.043 Rn.48779 G protein-coupled receptor AF281635 1 0.43 0.021 0.049 0.002 Rn.9264 zinc finger protein 22 (KOX 15) U24174 1 0.09 0.021 0.219 0.133 Rn.10089 cyclin-dependent kinase inhibitor 1A AW142947 1 -0.22 0.019 -0.085 -0.030 Rn.61563 ESTs BF396132 1 -0.26 0.014 -0.055 0.004 Rn.76362 echinoderm microtubula associated protein like in the composition of the	
AI102732 1 -0.31 0.035 -0.113 0.064 Rn.7539 ESTs AI176933 1 0.46 0.035 0.076 -0.048 Rn.23658 ajuba (Mm.) (DBSS) AF208288 1 -0.27 0.034 -0.127 0.043 Rn.48779 G protein-coupled receptor 26 AF281635 1 0.43 0.021 0.049 0.002 Rn.9264 zinc finger protein 22 (KOX 15) U24174 1 0.09 0.021 0.219 0.133 Rn.10089 cyclin-dependent kinase inhibitor 1A AW142947 1 -0.22 0.019 -0.085 -0.030 Rn.61563 ESTs BF396132 1 -0.26 0.014 -0.055 0.004 Rn.76362 echinoderm microtubule associated protein like in a second control of the composition o	1
AI176933 1 0.46 0.035 0.076 -0.048 Rn.23658 ajuba (Mm.) (DBSS) AF208288 1 -0.27 0.034 -0.127 0.043 Rn.48779 G protein-coupled receptor AF281635 1 0.43 0.021 0.049 0.002 Rn.9264 zinc finger protein 22 (KOX 15) U24174 1 0.09 0.021 0.219 0.133 Rn.10089 cyclin-dependent kinase inhibitor 1A AW142947 1 -0.22 0.019 -0.085 -0.030 Rn.61563 ESTs BF396132 1 -0.26 0.014 -0.055 0.004 Rn.76362 echinoderm microtubule associated protein like in the interval of the inter	
AF281635 1 0.43 0.021 0.049 0.002 Rn.9264 zinc finger protein 22 (KOX 15) U24174 1 0.09 0.021 0.219 0.133 Rn.10089 cyclin-dependent kinase inhibitor 1A AW142947 1 -0.22 0.019 -0.085 -0.030 Rn.61563 ESTs BF396132 1 -0.26 0.014 -0.055 0.004 Rn.76362 echinodern microtubul associated protein like sassociated protein like subscited protein like	
U24174	otor
AW142947 1 -0.22 0.019 -0.085 -0.030 Rn.61563 ESTs BF396132 1 -0.26 0.014 -0.055 0.004 Rn.76362 echinodern microtubul associated protein like in associated protein like in a subfactor receptor U57049 1 -0.17 0.013 -0.080 0.000 Rn.10494 reductase AW520754 1 -0.08 0.010 -0.124 0.021 Rn.15536 potassium channel, subfamily K, member 3	
BF396132	:
NM_012610 1 -0.08 0.014 -0.164 0.054 Rn.10980 nerve growth factor receptor	
NM_012610	
U57049 1 -0.17 0.013 -0.080 0.000 Rn.10494 methylenetetrahydrofola reductase AW520754 1 -0.08 0.010 -0.124 0.021 Rn.15536 potassium channel, subfamily K, member 3)
AW520754 1 -0.08 0.010 -0.124 0.021 Rn.15536 potassium channel, subfamily K, member 3	ıte
(Hs.) (DBSS)	
AI231846 1 -0.13 0.008 -0.059 0.032 Rn.27 ESTs	
BE116947 1 0.05 0.006 0.126 -0.078 Rn.8045 ESTs	
AW917933 1 -0.04 0.005 -0.124 0.039 Rn.28424 ESTs AW144517 1 -0.05 0.005 -0.097 -0.004 Rn.13780 ESTs	
AW144517 1 -0.05 0.005 -0.097 -0.004 Rn.13780 ESTs AW920818 1 0.03 0.005 0.177 -0.078 Rn.11702 macrophage activation	,
AB021980 1 -0.05 0.003 -0.057 0.054 Rn.32872 delta-6 fatty acid	٥
AF087454 1 -0.29 0.001 -0.004 0.033 Rn.30019 potassium voltage-gate	ı
channel, subfamily Q, member 3	
BE097309 1 0.41 0.000 0.001 0.004 Rn.46694 Peregrin (Bromodomain and PHD finger-contain protein 1) (Hs.) (DBSS_strong)	
AW919837 1 -0.05 0.000 0.010 0.042 Rn.23432 adrenergic, alpha-2A-, receptor (Hs.) (DBSS)	
NM_013197 1 0.03 -0.007 -0.259 -0.286 Rn.32517 aminolevulinic acid synthase 2	
BF396955 1 0.77 -0.050 -0.065 -0.228 Rn.41236 PC4035 cell-cycle-dependent 350K nuclea	
BF281149 1 1.34 -0.057 -0.042 -0.226 Rn.3137 Hypothetical protein (Hs.) (DBSS_v Hypothetical protein (KIAA0008 (Hs.) (DBSS_weak)	eak)
AI412011 2 3.38 0.279 0.082 0.005 Rn.3738 RIKEN cDNA 0610012G03; expressed sequence AI839730 (M (DBSS_weak)	
BF419406 2 -0.94 0.159 -0.168 -0.026 Rn.26560 ESTs	
NM_021682 2 -0.53 0.125 -0.234 -0.032 Rn.42884 kilon	
AF136583 2 0.66 0.115 0.174 -0.024 Rn.12100 serum-inducible kinase	
NM_020308 2 0.94 0.111 0.118 -0.025 Rn.28393 a disintegrin and metalloproteinase doma	
(ADAM) 15 (metargidi BE109152 2 1.60 0.103 0.064 0.011 Rn.19642 Red protein (RER protein (Mm.) (DBSS_strong)	
AI176739 2 0.41 0.083 0.205 0.005 Rn.22359 KIAA1002 protein (Hs. (DBSS_strong) (IDBSS_strong) (IDBSS_strong) (IDBSS_moderate))
AI228233 2 0.67 0.076 0.113 -0.017 Rn.25139 epsin 2 (Hs.) (DBSS)	
AF007549 2 0.55 0.075 0.136 0.026 Rn.10734 golgi SNAP receptor complex member 2	
AI232347 2 -2.15 0.070 -0.032 0.012 Rn.102 chromosome 14 open reading frame 114 (Hs.)
AW915996 2 -0.48 0.054 -0.114 0.094 Rn.19250 T00260 hypothetical protein KIAA0605 (Hs. (DBSS_strong))

TABLE 4-continued

	186 gei	nes identi	fied to be	necessary an	nd sufficient	to classify	the training set.
Probe	Iteration	Weight	Impact	Mean Logratio Positive Class	Mean Logratio Negative Class	Unigene ID	UniGene Description
AA819832	2	-0.40	0.054	-0.136	0.141	Rn.34433	period homolog 1 (Drosophila) (Hs.) (DBSS)
AW524724	2	-0.34	0.052	-0.156	-0.002	Rn.95059	ryanodine receptor type 1 (Mm.) (DBSS_strong)
BE103916	2	-0.72	0.046	-0.064	0.020	Rn.26832	ESTs
BF283302	2	0.56	0.046	0.081	-0.008	Rn.226	ESTs
X68878	2	-0.17	0.040	-0.244	-0.050	Rn.11022	synaptosomal-associated protein, 91 kDa
D00403	2	-0.44	0.039	-0.088	0.031	Rn.12300	Interleukin 1 alpha
AI145385	2	-0.79	0.035	-0.044	-0.025	Rn.3580	ESTs
AI317854	2	-0.22	0.032	-0.143	0.012	Rn.20362	ESTs
A I231432	2	0.58	0.030	0.051	-0.025	Rn.6983	hypermethylated in cancer 1 (Mm.) (DBSS_moderate)
AA 996961	2	-0.34	0.029	-0.088	0.071	Rn.12469	DNA-repair protein complementing XP-A cells (Hs.) (DBSS_moderate)
NM_012971	2	-0.26	0.025	-0.098	0.058	Rn.9884	potassium voltage gated channel, shaker related subfamily, member 4
BF397726	2	0.43	0.020	0.047	-0.076	Rn.18639	NF-E2-related factor 2 (Rn.) (DBSS_weak)
AW527217	2	-0.20	0.017	-0.088	-0.027	Rn.23378	ESTs
AA799789	2	0.25	0.016	0.065	-0.026	Rn.30163	ESTs
NM_013190	2	-0.59	0.015	-0.026	0.001	Rn.4212	Phosphofructokinase, liver, B-type
AI576621	2	0.16	0.013	0.082	0.027	Rn.24920	ESTs
AA943149	2	0.81	0.010	0.012	-0.002	Rn.7346	ALEX3 protein (Hs.) (DBSS_strong)
AW253895	2	-0.12	0.006	-0.055	0.011	Rn.3382	BRCA1 associated protein- 1 (ubiquitin carboxy- terminal hydrolase) (Hs.) (DBSS_strong)
BF283340	2	-0.09	0.005	0.057	0.028	Rn.20857	ESTs
AF073379	2	-0.11	0.005	-0.057 -0.046	0.028	Rn.10169	glutamate receptor, ionotropic, N-methyl-D- aspartate 3A
AA799981	2	-0.14	0.005	-0.034	0.032	Rn.6263	ESTs
AF237778	2	-0.18	0.003	-0.017	0.086	Rn.88349	calcium/calm odulin- dependent protein kinase II alpha subunit
AI175375	2	-0.14	0.003	-0.019	-0.025	Rn.24087	ESTs
A J130946	2	0.13	0.002	0.014	-0.096	Rn.2949	karyopherin (importin) alpha 2
AI012120	2	0.25	-0.004	-0.016	-0.149	Rn.17809	ESTs
AW252871	2	0.54	-0.078	-0.145	-0.370	Rn.12774	cell proliferation antigen Ki-67 (Mm.)
102.062	2	0.70	0.163	0.222	0.200	D 0010	(DBSS_moderate)
103863 Ј19614	3	0.70 2.55	0.163 0.161	0.233 0.063	0.208 -0.005	Rn.9918 Rn.11373	serine dehydratase lamina-associated polypeptide 1C
M19651	3	0.78	0.131	0.168	0.052	Rn.11306	Fos-like antigen 1
A I407719	3	-1.78	0.111	-0.063	0.161	Rn.20359	ubiquitin specific protease 2 (Hs.) (DBSS)
BF396629	3	2.54	0.111	0.044	-0.051	Rn.16544	patched homolog (Drosophila) (Hs.) (DBSS)
BF290678	3	2.25	0.109	0.049	-0.015	Rn.40449	heterogeneous nuclear ribonucleoprotein G (Mm.) (DBSS)
BE101099	3	-1.84	0.109	-0.059	-0.008	Rn.35019	parathyroid hormone regulated sequence (215 bp
AI070303	3	-1.13	0.098	-0.086	0.019	Rn.21284	pancreasin (Hs.) (DBSS_moderate)
AA925559	3	-1.06	0.078	-0.074	0.031	Rn.25196	RIKEN cDNA 2610027L16 [(Mm.) (DBSS_strong)
AB005549	3	0.58	0.056	0.097	-0.026	Rn.31803	three-PDZ containing protein similar to <i>C.</i> elegans PAR3 (partitioning defect)

TABLE 4-continued

	186 gei	nes identi	fied to be	necessary a	nd sufficient	to classify	the training set.
Probe	Iteration	Weight	Impact	Mean Logratio Positive Class	Mean Logratio Negative Class	Unigene ID	UniGene Description
AI717140	3	-0.59	0.043	-0.072	-0.001	Rn.22400	ESTs
AA858817	3	-0.23	0.040	-0.171	0.079	Rn.22047	T46271 hypothetical protein DKFZp564P1263.1
BF284897	3	0.54	0.035	0.064	0.027	Rn.18772	(Hs.) (DBSS_moderate) hypothetical protein FLJ10579 (Hs.)
							(DBSS_moderate)
AW914881	3	0.27	0.034	0.123	0.036	Rn.22383	ESTs
BE106459 BF283556	3	-0.21 -0.14	0.033 0.027	-0.157 -0.188	-0.037 0.019	Rn.20259 Rn.7829	ESTs Homo sapiens clone 23785
BF283330	3	-0.14	0.027	-0.188	0.019	KII. / 829	mRNA sequence (Hs.) (DBSS)
M63282	3	0.31	0.016	0.050	0.084	Rn.9664	Activating transcription factor 3
AW533663	3	0.08	0.014	0.174	0.124	Rn.41672	Proline oxidase, mitochondrial precursor
L19656	3	-0.92	0.013	-0.014	0.048	Rn.10552	(Mm.) (DBSS_strong) 5-hydroxytryptamine (serotonin) receptor 6
NM_012852	3	0.11	0.009	0.083	-0.008	Rn.34834	5-Hydroxytryptamine (serotonin) receptor ID
AA946230	3	-0.22	0.008	-0.039	-0.023	Rn.47222	ESTs
BF405135	3	-0.36	0.008	-0.022	0.018	Rn.51262	ESTs
AA818949	3	-0.14	0.007	-0.052	0.002	Rn.20419	DnaJ homolog subfamily B member 12 (Hs.) (DBSS_moderate)
X79860	3	-0.36	0.006	-0.017	0.066	Rn.65877	H1SHR mRNA
AW253907	3	-0.08	0.005	-0.064	0.066	Rn.98601	ESTs
X89603	3	0.05	0.004	0.091	-0.049	Rn.11325	metallothionein 3
AA858649	3	-0.50	-0.002	0.004	0.004	Rn.16864	chromosome 13 open reading frame 9 (Hs.) (DBSS_strong)
AW529588	3	0.61	-0.003	-0.005	-0.040	Rn.28180	ESTs
BF550800	3	0.16	-0.004	-0.023	-0.307	Rn.36317	ESTs
BE111296	3	0.18	-0.014	-0.079	-0.174	Rn.19339	ESTs
AI113104	3	1.77	-0.086	-0.048	-0.262	Rn.12343	protein regulator of cytokinesis 1 (Hs.) (DBSS_moderate)
U53706	4	-1.14	0.159	-0.139	-0.021	Rn.10288	mevalonate pyrophosphate decarboxylase
L36459	4	0.89	0.152	0.171	-0.036	Rn.10045	Interleukin 9 receptor
BF410042	4	4.02	0.151	0.038	-0.030	Rn.31227	cardiac lineage protein 1 (Mm.) (DBSS)
AW915655	4	-2.26	0.129	-0.057 -0.096	0.000	Rn.14962	ESTs
AA944518 NM_012939	4 4	-1.07 -0.19	0.102 0.079	-0.096 -0.408	0.019 -0.002	Rn.34351 Rn.1997	ESTs Cathepsin H
BF408867	4	-0.19	0.079	-0.408	0.002	Rn.35618	mitochondrial translational
		57	003		015		release factor 1-like (Hs.) (DBSS_moderate)
AW915454	4	-0.26	0.052	-0.204	-0.028	Rn.14822	ESTs
BE113132	4	-0.37	0.042	-0.112	0.124	Rn.22381	guanine nucleotide exchange factor for Rap1; M-Ras-regulated GEF (Hs.) (DBSS)
AW143273	4	0.72	0.040	0.056	-0.020	Rn.11888	Rec8p, a meiotic recombination and sister chromatid cohesion phosphoprotein of the rad21p family (Hs.)
ATTIO15107		0.70	0.020	0.055	0.000	D 10002	(DBSS)
AW915107	4	0.70	0.039	0.055	-0.023	Rn.19003	ESTs
BE110577	4 4	0.96 0.39	0.038 0.034	0.040	-0.008 -0.008	Rn.14584 Rn.13195	ESTs ATP-binding cassette, sub-
AW141985	4	0.39	0.034	0.088	-0.008	кп.13193	family C (CFTR/MRP), member 4
AW140530	4	-0.35	0.029	-0.083	0.005	Rn.7679	tumor susceptibility protein 101 (tsg101) gene (Mm.) (DBSS)

TABLE 4-continued

	186 gei	nes identi	fied to be	necessary a	nd sufficient	to classify	the training set.
Probe	Iteration	Weight	Impact	Mean Logratio Positive Class	Mean Logratio Negative Class	Unigene ID	UniGene Description
BF420720	4	-0.31	0.026	-0.083	0.030	Rn.23998	ESTs
AW144399	4	-0.78	0.025	-0.032	0.068	Rn.15255	hypothetical protein FLJ10652 (Hs.) (DBSS_moderate)
AI411605	4	-0.30	0.024	-0.079	-0.095	Rn.20056	ESTs
NM_019123 AW920802	4 4	0.38 0.50	0.021 0.019	0.055 0.037	-0.025 -0.021	Rn.88072 Rn.36609	sialyltransferase 7c ribosomal protein L5 (Hs.) (DBSS)
AI228598	4	-0.70	0.018	-0.026	0.036	Rn.11771	ESTs
AI175454	4	0.18	0.013	0.072	-0.002	Rn.17244	procollagen-proline, 2- oxoglutarate 4-dioxygenase (proline 4-hydroxylase), alpha polypeptide II (Hs.) (DBSS_strong)
AI009623	4	-0.08	0.011	-0.135	-0.073	Rn.13924	ESTs
AI235282	4	-0.20	0.011	-0.053	0.004	Rn.22436	Low-density lipoprotein receptor-related protein 1 precursor (Hs.) (DBSS_strong)
NM_012564	4	-0.06	0.009	-0.159	-0.100	Rn.1437	Group-specific component (vitamin D-binding protein)
BE095865	4	-0.35	0.009	-0.025	0.104	Rn.21852	calcium channel, voltage- dependent, alpha 1I subunit (Hs.) (DBSS)
AF291437	4	-0.40	0.009	-0.022	-0.058	Rn.39124	leucine rich repeat protein 3, neuronal
AF176351 AB027155	4	-0.26 0.15	0.009	-0.032 0.057	0.017	Rn.54003 Rn.44869	nuclear receptor coactivator 6 phosphodiesterase 10A
BE116569	4	0.34	0.008	0.024	-0.009	Rn.15835	zinc-finger protein AY163807 (Hs.) (DBSS_strong)
AA894210 AJ237852	4 4	0.05 -0.04	0.004 0.003	0.091 -0.058	0.082 0.065	Rn.85480 Rn.30023	ESTs sodium channel, voltage- gated, type1 1, alpha polypeptide
AJ305049	4	-1.09	0.002	-0.002	0.075	Rn.64632	interleukin 10 receptor, alpha
NM_017186	4	-0.03	0.002	-0.070	-0.015	Rn.30042	glial cells missing (Drosophila) homolog a
AA800004	4	0.04	0.001	0.024	-0.063	Rn.6269	Septin 4 (Peanut-like protein 2) (Brain protein H5) (Hs.) (DBSS_strong)
NM_012614 BF285985	4 4	0.05 -0.06	0.001 -0.001	0.012 0.016	0.040 0.074	Rn.9714 Rn.42366	Neuropeptide Y protein tyrosine phosphatase, receptor type, f polypeptide (PTPRF), interacting protein (liprin), alpha 4
AI412889	4	-0.08	-0.001	0.012	0.105	Rn.23659	monocyte to macrophage differentiation-associated 2 (Mm.) (DBSS)
AJ002556	4	-0.54	-0.003	0.006	0.050	Rn.37490	microtubule-associated protein 6
AI179459	4	0.12	-0.011	-0.094	-0.152	Rn.31366	Kell blood group (Mm.) (DBSS_moderate)
AI603128	4	0.15	-0.019	-0.127	-0.330	Rn.13094	Cyclin A2 (Cyclin A) (Mm.) (DBSS_strong)
BE111688 NM_012892	4 5	1.72 -0.70	-0.082 0.128	-0.048 -0.184	-0.343 -0.127	Rn.23351 Rn.37523	cyclin B2 (Hs.) (DBSS_strong) amiloride-sensitive cation
NM_012892 BE098463	5	2.30	0.128	-0.184 0.044	-0.127	Rn.3/523	channel 1 ESTs
C06844	5	-0.94	0.095	-0.101	0.075	Rn.7159	S49158 complement protein C1q beta chain precursor (Rn.) (DBSS_weak)

TABLE 4-continued

	186 gei	nes identi	fied to be	necessary ar	nd sufficient	to classify	the training set.
Probe	Iteration	Weight	Impact	Mean Logratio Positive Class	Mean Logratio Negative Class	Unigene ID	UniGene Description
AI170114	5	-0.42	0.078	-0.183	-0.112	Rn.91697	ESTs
AI105265	5	-1.53	0.073	-0.048	0.009	Rn.5911	hypothetical protein FLJ10315 (Hs.) (DBSS_strong)
BF394214	5	-0.79	0.071	-0.090	-0.014	Rn.58227	ESTs
AA946356	5	-1.08	0.063	-0.058	-0.017	Rn.1435	CGG triplet repeat binding protein 1 (Hs.) (DBSS)
AW919159	5	1.09	0.056	0.051	-0.022	Rn.41574	A38135 ADP- ribosylarginine hydrolase (Rn.) (DBSS_weak)
AI230884	5	1.61	0.053	0.033	-0.034	Rn.9797	Fibroblast growth factor receptor 1
BF406522	5	0.92	0.052	0.056	-0.019	Rn.3537	cerebellar degeneration- related protein 2, 62 kDa (Hs.) (DBSS)
NM_012848	5	0.14	0.048	0.350	0.110	Rn.54447	ferritin, heavy polypeptide 1
AW914090	5	-1.61	0.046	-0.029	0.002	Rn.973	60S acidic ribosomal protein P1 (Rn.) (DBSS_strong)
AW142828	5	-0.65	0.044	-0.068	-0.034	Rn.23877	ESTs
AI705731	5	-0.95	0.040	-0.042	0.058	Rn.24919	transcription factor MTSG1
NM_019126 U73503	5	-0.33 0.64	0.037	-0.112 0.057	0.140 -0.014	Rn.25723 Rn.10961	Carcinoembryonic antigen gene family (CGM3) calcium/calmodulin-
073303		0.04	0.037	0.037	-0.014	Kii.10501	dependent protein kinase (CaM kinase) II gamma
AF017437	5	0.55	0.036	0.066	-0.010	Rn.7409	integrin-associated protein
NM_021869	5	-0.42	0.035	-0.083	0.057	Rn.1993	syntaxin 7
AI144644	5	-0.34	0.030	-0.087	0.024	Rn.12319	ESTs
AA818377	5	0.79	0.029	0.037	-0.033	Rn.34063	hypothetical protein FLJ22419 (Hs.) (DBSS_weak)
AI171994	5	0.13	0.027	0.198	0.008	Rn.22380	ESTs
AA925167	5	-0.12	0.022	-0.180	0.106	Rn.8672	ESTs
BF398051	5	-0.38	0.020	-0.053	0.080	Rn.97322	ESTs
AW144075	5	0.48	0.019	0.040	-0.024	Rn.19790	ESTs
U26686	5	-0.09	0.015	-0.158	-0.045	Rn.10400	nitric oxide synthase 2
BF404426 U31866	5 5	-0.07 0.24	0.009 0.007	-0.128 0.029	-0.032 -0.037	Rn.63325 Rn.32307	ESTs Nclone10 mRNA
AW917475	5	-0.07	0.006	-0.087	0.055	Rn.16643	high-affinity immunoglobulin gamma Fc receptor I
AI408517	5	0.44	0.006	0.013	0.021	Rn.2773	protein phosphatase 1, regulatory (inhibitor) 5 subunit 14B
AF207605 AI178922	5 5	-0.34 -0.41	0.005 0.005	-0.015 -0.012	0.000 -0.023	Rn.42674 Rn.18670	tubulin tyrosine ligase leucine zipper and CTNNBIP1 domain containing (Hs.)
BF398403	5	0.41	0.005	0.011	-0.037	Rn.20421	(DBSS_moderate) mannosyl-oligosaccharide 1,3-1,6-alpha-mannosidase (EC 3.2.1.114) (Mm.) (DBSS_moderate)
M22923	5	0.05	0.004	0.091	-0.019	Rn.10922	membrane-spanning 4-domains, subfamily A, member 2
BE107747 BF281697	5 5	-0.05 0.57	0.004 0.004	-0.077 0.007	0.041 -0.024	Rn.29176 Rn.7770	ESTs potassium voltage-gated channel, Isk-related family, member 1-like (Hs.)
AB006461 AF100960	5 5	0.03 0.03	0.002 0.001	0.059 0.051	-0.009 -0.038	Rn.5653 Rn.8633	(DBSS) neurochondrin FAT tumor suppressor
							(Drosophila) homolog

TABLE 4-continued

186 genes identified to be necessary and sufficient to classify the training set.											
Probe	Iteration	Weight	Impact	Mean Logratio Positive Class	Mean Logratio Negative Class	Unigene ID	UniGene Description				
NM_017353	5	-0.21	-0.004	0.019	0.045	Rn.32261	tumor-associated protein 1				
AI231716	5	1.81	-0.007	-0.004	-0.138	Rn.24598	ESTs				
NM_012964	5	0.67	-0.024	-0.036	-0.298	Rn.92304	Hyaluronan mediated motility receptor (RHAMM)				
L06040	5	0.19	-0.035	-0.183	-0.306	Rn.11318	arachidonate 12- lipoxygenase				

[0156] The 186 genes of the necessary set listed in Table 4 correspond to 164 reward genes, of which 72 are induced on average across the nephrotoxicants. Additional genes not necessary for classification, but nonetheless differentially regulated by the nephrotoxicants relative to the negative class, were also considered.

Example 5

Using a Necessary Set to Generate New Signatures for Renal Tubule Injury

[0157] As shown above in Examples 1-3, a predictive signature for renal tubule injury comprising 35 genes may be derived using gene expression data from a microarray in the context of a chemogenomic database. Using the signature stripping method described above, four additional high performing predictive signatures for renal tubule injury may also be derived wherein each of the signatures is non-overlapping, i.e., comprises genes not used in any of the other signatures. Together, the union of the genes in these five signatures comprises a set of 186 genes that is necessary for deriving a predictive signature for renal tubule injury capable of classifying the training set above a selected threshold level of LOR=1.64.

[0158] This example demonstrates that additional signatures for renal tubule injury may be generated based on the necessary set of 186 genes. In addition, it is shown that at least four genes must be selected from the necessary set in order to generate a signature for renal tubule injury capable of performing above a selected threshold LOR of 4.00.

[0159] As listed in Table 4, for each gene from the necessary set of 186, an impact factor was calculated, corresponding to the product of the gene's weight and the

gene's expression mean logratio in the positive class (i.e., nephrotoxicants). Subsets of genes were chosen randomly from the necessary set of 186 so that the sum of the impacts of all genes in the subset accounted for 1, 2, 4, 8, 16, 32, or 64% of the total impact. Total impact was defined as the sum of the individual impacts of all 186 genes in the necessary set. This random subset selection procedure was repeated 20 times resulting in 140 gene subsets (i.e., 7 impact thresholds times 20 random choices).

[0160] Table 5 shows the average number of genes for each of these seven impact thresholds. This number increases regularly reaching an average of 116 genes for those subsets that account for 64% of the total impact. Each of these random subsets was used as input to compute a renal tubule injury signature using the SPLP algorithm as described in Example 3 above. A training LOR and a 10-fold cross-validated test LOR were calculated for each signature. Table 5 lists average LOR values for the signatures generated in each of the seven percent of total impact thresholds.

[0161] Based on the results tabulated in Table 5 it may be concluded that signatures for renal tubule injury capable of performing with an average training LOR of 4.30 may be generated starting with random subsets having an average of 4.4 genes that together have only 2% of the total impact of the necessary set. Similarly signatures capable of performing with an average test LOR of 4.41 may be derived from random subsets of the necessary set having an average of 9.15 genes with only 4% of the total impact. Significantly, the average training LOR never drops below 4.00 when a random set of genes having at least 4% impact are selected. As shown in Table 5, comparably higher performing signatures are derived from the necessary set when the random subsets have a percent impact of 8% or higher.

TABLE 5

RTI signatures generated based on randomly selecting necessary set genes with minimal percentage impact														
	# input genes Signature Length LOR (training) I													
percent impact*	avg	min	max	avg	min	max	avg	stdev	avg	stdev				
1	2.85	1	5	2.8	1	5	3.42	1.61	3.01	1.34				
2	4.4	1	9	4.3	1	8	4.30	1.61	3.20	1.00				
4	9.15	3	17	8.05	3	13	6.82	2.34	4.41	2.43				
8	17.3	8	27	12.8	8	18	8.54	0.61	5.91	1.99				
16	33.4	22	42	19.2	14	25	8.68	0.00	7.85	2.01				

TABLE 5-continued

RTI signatures generated based on randomly selecting necessary set genes with minimal percentage impact

	# in	put gen	.es	Sign	ature L	ength	LOR (training)	LOR	(test)
percent impact*	avg	min	max	avg	min	max	avg	stdev	avg	stdev
32 64	61.6 116	49 100	76 134	26.5 30.7	22 28	30 36	8.68 8.68	0.00 0.00	7.35 7.07	2.03 1.50

^{*}average of 20 lists chosen from the necessary set

[0162] Table 6 shows the parameters for 20 signatures generated from random subsets of genes with 2% of the total impact of the 186 gene necessary set. Tables 7 (subset 8) and 8 (subset 14) illustrate two specific 5 gene signatures (including values for gene weights and bias) for predicting renal tubule injury onset that perform with a training LOR of 4.00 and 7.3, respectively.

TABLE 6

RTI signatures generated based on random selections of necessary set genes with 2% impact

Subset #	# Input Genes	Signature Length	Training LOR	Test LOR
14	5	5	7.3	5.0
9	7	7	6.8	3.4
15	5	5	6.2	4.1
7	6	6	6.0	3.2
18	5	5	5.8	3.7
3	4	4	5.5	4.0
10	9	8	5.0	2.8
2	4	3	4.7	1.7
13	3	3	4.5	3.2
19	6	6	4.4	2.6
8	5	5	4.0	2.8
11	5	5	3.8	4.5
4	4	4	3.8	4.0
12	4	4	3.8	5.1
20	4	4	3.2	2.7
5	3	3	2.8	2.6
1	4	4	2.6	2.4
17	3	3	2.2	2.4
6	1	1	2.1	1.6
16	1	1	1.7	2.3

[0163]

TABLE 7

Subset	8	
BF283302	15.5	
AW920818	5.88	
AW141985	5.48	
BF403410	4.28	
AA858649	-2.3	
Bias	1.13	

[0164]

TABLE 8

	Subset 14
AI176933	43.1
U08257	33.7
BE116947	18.4
AI408517	12.7
AA819832	-2.9
Bias	8.49

[0165] Similarly Table 9 shows the parameters for 20 signatures generated from random subsets of genes with 4% of the total impact of the 186 gene necessary set. Tables 10 (subset 18) and 11 (subset 5) illustrate specific 9 and 13 gene signatures for predicting renal tubule injury onset that perform with a test LOR of 4.1 and 10.2, respectively.

TABLE 9

Subset #	# Input Genes	Signature Length	Training LOR	Test LOR
5	13	13	8.7	10.2
2	14	11	8.7	8.9
7	11	10	8.7	8.9
9	17	11	8.7	6.2
20	11	9	8.7	5.3
10	14	12	8.7	4.7
11	13	12	8.7	4.6
14	7	6	8.7	4.5
12	9	8	8.7	4.3
18	9	9	8.7	4.1
15	11	9	8.7	3.8
3	6	6	6.2	3.3
19	7	6	6.2	3.2
13	6	6	4.7	3.1
8	11	9	6.8	2.7
4	5	5	4.3	2.7
17	5	5	3.7	2.1
1	7	7	3.7	2.1
6	4	4	3.4	2.0
16	3	3	1.9	1.5

[0166]

TABLE 10

Subset	18
AW143273	55.95
AI599126	29.8
AI705731	19.05
BF406522	16.71

TABLE 10-continued

Sub	set 18	
AB027155	-4.12	
AW253895	-13.53	
AA819832	-14.81	
X68878	-17.57	
AW140530	-19.85	
Bias	8.96	

[0167]

TABLE 11

Subset 5	i	
AW144075	4.82	
AI113104	4.58	
AI171994	4.25	
AW920818	3.39	
BF281697	3.11	
AI012120	1.76	
BE110577	1.08	
NM 012964	0.87	
AI227912	0.74	
AW144399	-0.2	
AI232347	-2.9	
AA944518	-6.4	
AW914090	-6.6	
Bias	0.68	

[0168] The results tabulated in Table 5 may also be illustrated graphically. As shown in FIG. 2, which plots training LOR and test LOR versus signature length, a signature performing with an average training LOR of 4.00 may be achieved by randomly selecting on average 4 genes from the necessary set. Similarly, an average test LOR of 4.00 may be achieved by randomly selecting on average 7 genes from the necessary set.

Example 6

Functional Characterization of the Necessary Set of Genes for Renal Tubule Injury by Random Supplementation of a Fully Depleted Set

[0169] This example illustrates how the set of 186 genes necessary for classifying renal tubule injury may be functionally characterized by randomly supplementing and thereby restoring the ability of a depleted gene set to generate RTI signatures capable of performing on average above a threshold LOR. In addition to demonstrating the power of the 186 information rich genes in the RTI necessary set, this example illustrates a system for describing any necessary set of genes in terms of its performance parameters.

[0170] As described in Example 4, a necessary set of 186 genes (see Table 4) for the RTI classification question was generated via the stripping method. In the process, a corresponding fully depleted set of 7292 genes (i.e., the full dataset of 7478 genes minus 186 genes) was also generated. The fully depleted set of 7292 genes was not able to generate an RTI signature capable of performing with a LOR greater than or equal to 1.28 (based on cross-validation using 40 random 80:20 training:test splits).

[0171] A further 186 genes were randomly removed from the fully depleted set. Then a randomly selected set including 10, 20, 40 or 80% of the genes from either: (a) the necessary set; or (b) the set of 186 randomly removed from the fully depleted set; is added back to the depleted set minus 186. The resulting "supplemented" depleted set was then used to generate an RTI signature, and the performance of this signature is cross-validated using 3 random 60:40 training:test splits. This process was repeated 20 times for each of the different percentage supplementations of genes from the necessary set and the random 186 genes removed from the original depleted set. Twenty cross-validated RTI signatures were obtained for each of the various percentage supplementations of the depleted set. Average LOR values were calculated based on the 20 signatures generated for each percentage supplementation.

[0172] Results

[0173] As shown in Table 12, supplementing the fully depleted set (minus random 186) with as few as 10% of the randomly chosen genes from the necessary set results in significantly improved performance for classifying RTI. The random 10% of genes selected from the depleted 186 yielded signatures performing with an avg. LOR=1.4. In contrast, supplementing the depleted set (minus random 186) with 10% from the necessary set yields RTI signatures performing with an avg. LOR=4.5 (based on 3-fold cross-validation using random 60:40 splits).

TABLE 12

Supplementation	Supplementation with random genes from necessary or depleted sets			
%	Necessary Set Avg. LOR	Depleted Set Avg. LOR		
10 20 40 80	4.51 4.93 4.73 4.10	1.43 2.32 2.63 3.28		

[0174] Although increasing the percentage of random "depleted" set genes used to supplement resulted in an increase in average performance, even at 80%, the average LOR remained below 4.00, while supplementation with the random 80% "necessary" set genes yielded an average LOR above 4.00.

[0175] These results demonstrate how supplementation with a percentage of randomly selected genes from the RTI necessary set of 186 "revives" the performance of a fully depleted set for generating classifiers. Thus, the RTI necessary set of genes may be functionally characterized as the set of genes for which a randomly selected 10% will supplement a set of genes fully depleted for RTI classification (i.e., not capable of producing RTI signatures with avg. LOR>~1.4), such that the resulting "revived" gene set generates RTI signatures with an average LOR greater than or equal to 4.00.

Example 7

Functional Characterization of the RTI Necessary Set by Random Supplementation with Rigorous Signature Cross-Validation

[0176] In a further exemplification of the method of Example 6, a randomly selected set including 1, 2, 5, 10, 20,

40, 80, 90, or 99% of the genes from either: (a) the necessary set; or (b) the set of 186 randomly removed from the fully depleted set; was added back to the depleted set minus 186. The resulting "supplemented" depleted set was then used to generate an RTI signature, and the performance of this signature was cross-validated using 40 random 80:20 training:test splits. This process was repeated 100 times for each of the different percentage supplementations of genes from (a) the necessary set, and (b) the random 186 genes removed from the original depleted set. Twenty cross-validated RTI signatures were obtained for each of the various percentage supplementations of the depleted set. Average LOR values were calculated based on the 20 signatures generated for each percentage supplementation.

[0177] Results

[0178] Based on cross-validation using 40 random 80:20 training:test splits, the fully depleted set of 7292 genes was not able to generate an RTI signature capable of performing with a LOR greater than or equal to 1.28. As shown in Table 13, supplementing the fully depleted set (minus random 186) with as few as 5% of the randomly chosen genes from the necessary set results in substantially improved performance for classifying RTI (avg. LOR ~2.2). In contrast, the random 5% of genes selected from the depleted 186 yielded signatures performing with an avg. LOR ~1.3. Significantly, increasing the percentage of random "depleted" set genes used to supplement did not result in an increase in average performance—even at 99%, the average LOR remained at ~1.3, while supplementation with the random 99% "necessary" set genes yielded an average LOR of ~4.3.

TABLE 13

upplementation with r	andom genes from nec	pplementation with random genes from necessary or depleted set		
% Supplementation	Necessary Set Avg LOR	Random Set Avg LOR		
1	1.44	1.31		
2	1.72	1.31		
5	2.19	1.31		
10	2.68	1.31		
20	3.38	1.30		
40	4.00	1.30		
80	4.39	1.28		
90	4.32	1.28		
99	4.32	1.28		

[0179] These results further demonstrate how supplementation with even a small percentage of randomly selected genes from the RTI necessary set "revives" the performance of a fully depleted set for generating classifiers. It also demonstrates that more rigorous cross-validation (40-fold random 80:20 training:test splits) provides a more consistent average performance of the signatures generated by the random supplementations from depleted set. Thus, the RTI necessary set of genes may be functionally characterized as the set of genes for which a randomly selected 5% will supplement a set of genes fully depleted for RTI classification (i.e., not capable of producing RTI signatures with avg. LOR>~1.3), such that the resulting "revived" gene set generates RTI signatures with an average LOR of greater than or equal to about 2.00. Further, a random supplementation of at least 40% of the necessary set genes will produce a revived gene set capable of generating RTI signatures with an average LOR greater than or equal to about 4.00.

Example 8

Construction and Use of a DNA Array for Predicting Renal Tubule Injury

[0180] The necessary subset of 186 genes identified to be necessary and sufficient to classify the renal tubule injury training set listed in Table 4 may be used as the basis for a DNA array diagnostic device for predicting renal tubule injury. The device may be used in a therapeutic monitoring context, such as for monitoring the response of an individual to a compound that is suspected of possibly causing renal tubule injury (or related nephrotoxic side effects). Alternatively, smaller sufficient subsets of genes the necessary set, which may be selected according to the methods of Examples 4 and 5 described above, may be used as the basis for a DNA array.

[0181] The probe sequences used to represent the 186 (or fewer) genes on the array may be the same ones used on the Amersham CodeLink™ RU1 platform DNA array used to derive the renal tubule injury signature as described in Examples 1-3. The 186 probes are pre-synthesized in a standard oligonucleotide synthesizer and purified according to standard techniques. The pre-synthesized probes are then deposited onto treated glass slides according to standard methods for array spotting. For example, large numbers of slides, each containing the set of 186 probes, are prepared simultaneously using a robotic pen spotting device as described in U.S. Pat. No. 5,807,522. Alternatively, the 186 probes may be synthesized in situ one or more glass slides from nucleoside precursors according to standard methods well known in the art such as ink-jet deposition or photoactivated synthesis.

[0182] The DNA probe arrays made according to this method are then each hybridized with a fluorescently labeled nucleic acid sample. The nucleic acid may be derived from mRNA obtained from a biological fluid (e.g., blood) or a tissue sample from a compound treated individual. Any of the well-known methods for preparing labeled samples for DNA probe array hybridization may be used. The fluorescence intensity data from hybridization of the sample to the DNA array of 186 (or fewer) genes of the necessary set is used to calculate expression log ratios for each of the genes. Depending on the specific gene signature selected for use in predicting renal tubule injury (e.g., the genes in iteration 1 of Table 4), the scalar product for that signature is calculated (i.e., sum of the products of expression log₁₀ ratio and weight for each gene less the bias). If the scalar product is greater than zero then the sample is classified as positive (i.e., onset of renal tubule injury is predicted).

[0183] All publications and patent applications cited in this specification are herein incorporated by reference as if each individual publication or patent application were specifically and individually indicated to be incorporated by reference.

[0184] Although the foregoing invention has been described in some detail by way of illustration and example for clarity and understanding, it will be readily apparent to one of ordinary skill in the art in light of the teachings of this invention that certain changes and modifications may be made thereto without departing from the spirit and scope of the appended claims.

What is claimed:

- 1. A reagent set for testing whether renal tubule injury will occur in a test subject comprising a plurality of polynucleotides or polypeptides representing a plurality of genes selected from Table 4.
- 2. The reagent set of claim 1, wherein the plurality of genes is the set of genes in any one of iterations 1 through 5 in Table 4.
- 3. The reagent set of claim 1, wherein the plurality of genes are selected from a linear classifier capable of classifying renal tubule injury with a training log odds ratio of greater than or equal to 4.35.
- **4**. The reagent set of claim 1, wherein the plurality of genes includes at least 4 genes having at least 2% of the total impact of all of the genes in Table 4.
- 5. The reagent set of claim 1, wherein the plurality of genes includes at least 8 genes having at least 4% of the total impact of the genes in Table 4.
- **6**. The reagent set of claim 1, wherein the reagents are polynucleotide probes capable of hybridizing to the plurality of genes selected from Table 4.
- 7. The reagent set of claim 6, wherein the polynucleotide probes are primers for amplification of the plurality of genes.
- **8**. The reagent set of claim 6, wherein the polynucleotide probes are immobilized on one or more solid surfaces.
- **9**. The reagent set of claim 1, wherein the reagents are polypeptides that bind to a plurality of proteins encoded by the plurality of genes selected from Table 4.
- 10. The reagent set of claim 9, wherein the proteins are secreted proteins.
- 11. An apparatus for predicting whether renal tubule injury will occur in a test subject comprising a reagent set according to claim 1.
- 12. The apparatus of claim 11, wherein the reagents are polynucleotides.

- 13. The apparatus of claim 11, wherein the reagents are polypeptides.
- 14. A set of genes useful for testing whether a compound will induce renal tubule injury comprising a random selection of at least about 10% of the genes from Table 4, wherein the addition of said randomly selected genes to a fully depleted gene set for the renal tubule injury classification question increases the average logodds ratio of the linear classifiers generated by the depleted set to at least about 2.5.
- 15. The set of claim 14, wherein the randomly selected percentage of genes from the necessary set is at least 20% and the average logodds ratio is increased to at least about 3.3
- 16. The set of claim 14, wherein the randomly selected percentage of genes from the necessary set is at least 40% and the average logodds ratio is increased to at least about 4.0.
- 17. A reagent set for classifying renal tubule injury comprising a set of polynucleotides or polypeptides representing a plurality of genes selected from Table 4, wherein the addition of a random selection of at least 10% of said plurality of genes to the fully depleted set for the renal tubule injury classification question increases the average logodds ratio of the linear classifiers generated by the depleted set by at least 2-fold.
- 18. The reagent set of claim 17, wherein the random selection is of at least 40% of said plurality of genes and the average logodds ratio of the linear classifiers generated by the depleted set by at least 3-fold.
- 19. An apparatus comprising a set of polynucleotides capable of specifically binding to the reagent set of claim 17.

* * * * *