

(19) **DANMARK**

(10) **DK/EP 3149032 T3**



(12)

Oversættelse af europæisk patentskrift

Patent- og
Varemærkestyrelsen

-
- (51) Int.Cl.: **C 07 K 16/00 (2006.01)** **C 07 K 14/39 (2006.01)** **C 07 K 16/14 (2006.01)**
C 07 K 16/46 (2006.01)
- (45) Oversættelsen bekendtgjort den: **2021-08-30**
- (80) Dato for Den Europæiske Patentmyndigheds bekendtgørelse om meddelelse af patentet: **2021-07-07**
- (86) Europæisk ansøgning nr.: **15726923.4**
- (86) Europæisk indleveringsdag: **2015-05-28**
- (87) Den europæiske ansøgnings publiceringsdag: **2017-04-05**
- (86) International ansøgning nr.: **EP2015061819**
- (87) Internationalt publikationsnr.: **WO2015181282**
- (30) Prioritet: **2014-05-29 GB 201409558**
- (84) Designerede stater: **AL AT BE BG CH CY CZ DE DK EE ES FI FR GB GR HR HU IE IS IT LI LT LU LV MC MK MT NL NO PL PT RO RS SE SI SK SM TR**
- (73) Patenthaver: **UCB Biopharma SRL, Allée de la Recherche 60, 1070 Brussels, Belgien**
- (72) Opfinder: **FINNEY, Helene Margaret, IPD, UCB Celltech, 208 Bath Road, Slough, Berkshire SL1 3WE, Storbritannien**
RAPECKI, Stephen Edward, IPD, UCB Celltech, 208 Bath Road, Slough, Berkshire SL1 3WE, Storbritannien
WRIGHT, Michael John, IPD, UCB Celltech, 208 Bath Road, Slough, Berkshire SL1 3WE, Storbritannien
- (74) Fuldmægtig i Danmark: **Zacco Denmark A/S, Arne Jacobsens Allé 15, 2300 København S, Danmark**
- (54) Benævnelse: **Nyt bispecifikt format, der er egnet til anvendelse i high-through-put-screening**
- (56) Fremdragne publikationer:
WO-A1-93/11162
WO-A1-2004/081051
WO-A1-2014/001326
D. M. GOLDENBERG ET AL: "Multifunctional Antibodies by the Dock-and-Lock Method for Improved Cancer Imaging and Therapy by Pretargeting", THE JOURNAL OF NUCLEAR MEDICINE, vol. 49, no. 1, 1 January 2008 (2008-01-01), pages 158-163, XP055028152, ISSN: 0161-5505, DOI: 10.2967/jnumed.107.046185
ROLAND KONTERMANN: "Dual targeting strategies with bispecific antibodies", MABS, LANDES BIOSCIENCE, US, vol. 4, no. 2, 1 March 2012 (2012-03-01), pages 182-197, XP009160769, ISSN: 1942-0870, DOI: 10.4161/MABS.4.2.19000
None

Description**Field of Invention**

[0001] The present disclosure relates to a method, in particular an *in vitro/ex vivo* method, of detecting synergistic biological function in a heterodimerically-tethered bispecific protein complex, libraries/multiplexes of the bispecific protein complexes, and kits and compositions thereof. The disclosure further relates to said novel bispecific protein complexes and use of the same in therapy, research and experimental purposes (in particular in assays looking for synergistic biological function). The present disclosure also extends to methods of preparing said bispecific complexes.

Background of Invention

[0002] Biological mechanisms *in vivo* are extremely complicated cascades of signals, which are difficult to deconvolute and understand. An example of such signalling is that required to activate T-cells, see Figure 1, from www.cellsignal.com. Activation of T cells requires at least two signals.

[0003] The recognition of the antigen by the T cell receptor is considered the first signal and the second signal arises from co-stimulation which results from the ligation of additional surface molecules on the T cell with additional molecules on an antigen presenting cell.

[0004] Thus T cell activation can be used to illustrate that the modulation of biological functions can require multiple signals. Other biological processes are equally complicated or more complicated. Whilst *in vitro* screening based on cells has and can assist with gaining insights into *in vivo* mechanisms the problem still arises of how to identify appropriate ligand pairs which modulate the biological function.

[0005] Bispecific antibodies are widely expected to play a major role in the next generation of biotherapeutics (D. Holmes, Nature Rev Drug Disc Nov 2011:10; 798). They have the potential to deliver superior, long term, broad efficacy in a greater proportion of patients.

[0006] This can be achieved by either co-engaging different antigens simultaneously within a common disease pathway, thereby reducing redundancy; or by targeting antigens from independent pathways to provide an additive or synergistic effect.

[0007] Bispecific antibodies facilitate access to novel biology such as:

- 1) cross-linking receptors on a cell,
- 2) inducing cell mediated effects,
- 3) localizing a cytokine to a cell to regulate signaling or locally block cytokine function,
- 4) engaging multiple epitopes simultaneously to generate "new activity", increase function or specificity, which may not be exhibited by a single monoclonal antibody or indeed mixtures of un-linked antibodies ('poly-monoclonals').

[0008] Current strategies to engage dual targets are largely based on rational design of known mechanisms and include: cross-linking inhibitory receptors, co-engagement/clustering of receptors, blocking multiple stimulatory pathways, selective engagement of inhibitory receptors and blocking distinct pathways such as co-stimulation & cytokine signaling.

[0009] However, the current state of the art in relation to known mechanisms and targets is a limiting factor to progress in this area.

[0010] Whilst bispecific antibodies have enormous potential as biological therapeutics they also present an increased set of challenges within discovery and development compared to monoclonal antibodies. Two key areas of difficulty are, 1) the development of a successful bispecific antibody format, and 2) selecting the pairs of targets to which the bispecific antibody will crosslink or co-engage.

[0011] Many promising bispecific antibody formats have now been developed that could potentially work as successful therapeutics including DVD-Ig (Abbvie), DuoBodies (Genmab), Knobs-in-Holes (Genentech), Common light chain (Merus). However, in each of these cases these formats are not ideally suited to high throughput target-dual-antigen discovery screening to enable the discovery of novel antigen pairs for crosslinking with bispecific antibodies.

[0012] Typically for a single bispecific antibody construct at least two variable regions need to be sub-cloned from the original source of discovery vectors (e.g. phage display, hybridoma or single B-cell cloning) into appropriate bispecific expression vectors, each arm of the bispecific has to be expressed and the resulting bispecific antibody purified. This cloning and subsequent expression effort quickly becomes a significant practical bottleneck if large numbers of pairs of variable regions are to be combined in an attempt to screen for the most efficacious combination of discovered variable regions or to discover novel antigen pairs.

[0013] For example, if 50 unique antibodies are discovered against a panel of 50 cell surface targets, then a total of 2500 bispecific antibodies could potentially be generated (envisaged as an X-by-Y grid). With the bispecific antibody

formats described above this would require at least 100 individual cloning reactions (50-X and 50-Y) followed by 2500 antibody expression experiments. Increasing the number of starting monoclonal antibodies to 100 would increase the minimal number of cloning reactions to 200 (100-X and 100-Y) and the expression number to 10,000.

[0014] Generally the root cause of this 'expression bottleneck' is the fact that the formats described above require both protein chain 'halves' of the final bispecific construct to be expressed simultaneously within a single expression experiment in the same cell. Therefore, for many formats, to produce 2500 bispecific antibodies, 2500 expression experiments are required.

[0015] The 'expression bottleneck' is further exacerbated if the bispecific antibody format is monocistronic (i.e. cloned and expressed as a single chain protein), for example single chain diabodies, as the number of cloning experiments would be 2500 and 10,000 respectively for the numbers given above.

[0016] Furthermore after expression, extensive purification may be required to isolate the desired construct.

[0017] Some bispecific approaches employ a common light chain in the bispecific constructs in order to reduce the amount of cloning, although this doesn't reduce the number of expression experiments. Furthermore, using a common chain, such as a common light chain, makes the challenge of antibody discovery harder as it is more difficult to find the starting antibody variable domains as the antibody needs to bind its antigen with a high enough affinity through one chain, such as the heavy chain, alone.

[0018] Accordingly the use of current bispecific formats in large scale and high throughput screening to identify novel antigen pairs is impractical and has led to the continued use of solely hypothesis driven approaches to bispecific antigen targeting.

[0019] We propose that rather than designing and testing a limited selection of bispecific antibodies that engage given epitopes on two known targets, the true potential of exploiting access to novel biology with bispecific antibodies can only be achieved through a broad functional screening effort with a large, diverse combinatorial panel of bispecific antibodies or protein ligands. To facilitate this screening a format and a method is required that enables the generation of large numbers of diverse bispecific proteins which can be readily constructed and screened for functional effects in a variety of functional screens. This approach allows for the serendipitous identification of synergistic pairs.

[0020] Thus it would be useful to generate and screen a large number of bispecific protein complexes present as combinations of various antigen specificities. In particular, it would be useful to be able to generate and screen a large number of different bispecific antibody complexes in a quick and efficient manner. There are a range of existing methods for manufacturing bispecific antibodies as already described above. However, each of these methods has its disadvantages, as do alternative methods as further described in more detail below.

[0021] The problem of how to efficiently identify targets for bispecific and multispecific constructs has not been adequately addressed in the art. For example WO2014/001326 employs chemical conjugation of a protein to a DNA fragment, wherein the DNA fragment hybridises to a complementary DNA sequence that links two such proteins together for generating tailor-made patient-specific multispecific molecules comprising at least two targeting entities.

[0022] There are number of difficulties associated with this approach if it were to be applied to identifying new bispecific combinations, for example conjugation of the protein to the DNA can result in damage to the activity and/or structure of the protein. In particular protein-DNA hybrids are not naturally occurring thus there is a potential for interference. In addition the chemical conjugation required to join the protein and DNA adds complexity, time and expense to the process.

[0023] Coupling and conjugation techniques exist for generating antibody drug conjugates and *in vivo* targeting technologies. Traditional chemical cross-linking is labour intensive as the relevant species may need to be purified from homodimers and other undesirable byproducts. In addition, the chemical modification steps can alter the integrity of the proteins, thus leading to poor stability or altered biological function. As a result, the production of bispecific antibodies by chemical cross-linking is often inefficient and can also lead to a loss of antibody activity.

[0024] Another method of manufacturing bispecific antibodies is by cell-fusion (e.g. hybrid hybridomas), wherein the engineered cells express two heavy and two light antibody chains that assemble randomly. Since there are 4 possible variants to choose from, this results in the generation of 10 possible bispecific antibody combinations, of which only some (in many cases, only one) combinations would be desired. Hence, generating bispecific antibodies by cell-fusion results in low production yields and also requires an additional purification step in order to isolate the desired bispecific antibodies from the other bispecific antibodies produced. These disadvantages increase manufacturing time and costs.

[0025] Recombinant DNA techniques have also been employed for generating bispecific antibodies. For example, recombinant DNA techniques have also been used to generate 'knob into hole' bispecific antibodies. The 'knob into hole' technique involves engineering sterically complementary mutations in multimerization domains at the CH3 domain interface (see e.g., Ridgway et al., Protein Eng. 9:617-621 (1996); Merchant et al., Nat. Biotechnol. 16(7): 677-81 (1998); see also U.S. Pat. Nos. 5,731,168 and 7,183,076). One constraint of this strategy is that the light chains of the two parent antibodies have to be identical to prevent mispairing and formation of undesired and/or inactive molecules when expressed in the same cell. Each bispecific (heavy and light chains thereof) must be expressed in a single cell and the protein product generally contains about 20% of homodimer, which is subsequently removed by purification.

[0026] Other approaches are based on the natural exchange of chains in full-length IgG4 molecules (Genmab Dubody).

However, this approach also has difficulties because it does not allow a construct to be prepared without an Fc region. As the Fc region can contribute to biological activity it may be difficult to establish if an activity observed is based on the combination of variable regions, the Fc or both in bispecific molecules comprising an Fc. Furthermore, the exchange is a dynamic process and this may lead to difficulties in relation to what the entity tested actually is.

[0027] Finally, a leucine zippers' approach has been described in WO 93/11162. A bispecific antibody comprising an anti-Tac Fab linked to the jun leucine zipper (anti-Tac-Jun) and an anti-CD3 linked to the fos leucine zipper (anti-CD3-Fos) was created through the pairwise association of fos and jun, leading to the formation of a bispecific F(ab'-zipper)₂ dimer. Such approach, however, suffers the drawback that along heterodimers, monospecific homodimers (jun-jun or fos-fos) will also be formed.

[0028] Thus there is a need for new methods of generating bispecific protein complexes to enable the more efficient and higher throughput screening of bispecific antibodies. In particular, there is a need for a format and a method wherein a selection of any two antibodies or antibody fragments from a pool of available antibodies or antibody fragments can be readily combined to efficiently produce a multiplex of different bispecific antibodies, whilst, for example avoiding or minimising the formation of homodimers. Assembling different bispecific antibodies efficiently is particularly important when screening for synergistic biological function for new combinations of antigen specificities, in particular where heterodimers are essential for discovering that function.

Summary of Invention

[0029] In one aspect there is provided a new bispecific format particularly suitable for use in screening because all of the components can be expressed from a cell as individual units, essentially without aggregation and the units can be assembled simply by mixing without employing conjugation or coupling chemistry and with minimal homodimerisation.

[0030] Thus there is provided a bispecific protein complex having the formula A-X:Y-B wherein:

A-X is a first fusion protein;

Y-B is a second fusion protein;

X:Y is a heterodimeric-tether;

wherein X is independently selected from a scFv, or a peptide, with the proviso that when X is a peptide Y is a scFv and when X is scFv then Y is a peptide;

wherein X or Y is a scFv specific to the peptide GCN4 represented by SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1; wherein the scFv is designated 52SR4 and is represented by SEQ ID NO:3 or amino acids 1 to 243 of SEQ ID NO:3 and wherein X or Y is a peptide GCN4 (SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1).

[0031] Hence, in the present invention the variable X or Y is the 52SR4 scFv (SEQ ID NO:3 or amino acids 1 to 243 of SEQ ID NO:3) and the other variable of X and Y is the peptide GCN4 (SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1).

[0032] Thus the bispecific format comprises two Fab arms with different specificities linked, for example via their C-terminal, by a scFv and peptide binding interaction. This type of arrangement is ideal for use in screening because there is no difficulty expressing the unit AX or the unit B-Y. The Fab/Fab' fragment is very stable and is not susceptible to inappropriate dimerization. Thus the amount of purification required after expression of each unit (A-X or B-Y) is minimal or in fact, unnecessary. The bispecific complex can be formed in a 1:1 molar ratio by simply admixing the relevant units i.e. without recourse to conjugation and coupling chemistry. The constant regions in the Fab/Fab' fragment drive dimerization of the Fab/Fab' components and the binding partners X and Y drive the equilibrium further in favour of forming the requisite heterodimer bispecific complex. Again little or no purification is required after formation of the complex after heterodimerisation. Thus large number of A-X and B-Y can be readily prepared and combined.

[0033] The Fab/Fab' entities in the complex mean the binding domains are held in biologically relevant orientations which mimic classic antibody geometry and this may contribute to the success of translating the pairs of variable regions identified by the screening method described herein below into other bispecific therapeutic formats which retain activity. The ability to prepare and screen a bispecific complex lacking the Fc fragment CH₂-CH₃ also ensures that the biological activity observed is in fact due solely to the variable region pairs in the complex. The simplicity of the bispecific complex of the invention and the methods of preparing it are a huge advantage in the context of facilitating high-through-put screening of variable domain pairs to find new target antigen combinations and also to optimise variable region sequences for a given combination.

[0034] In one embodiment A is a Fab fragment. In one embodiment B is a Fab fragment. In one embodiment A and B are both a Fab fragment (also referred to herein as a Fab-Kd-Fab).

[0035] In one embodiment X is fused to the C-terminal of a heavy chain or light chain in the Fab or Fab' fragment, in particular the C-terminal of the heavy chain.

[0036] In one embodiment Y is fused to the C-terminal of a heavy chain or light chain in the Fab or Fab' fragment, in

particular the C-terminal of the heavy chain.

[0037] X is independently selected from the 52SR4 scFv, and the GCN4 peptide, with the proviso that when X is the peptide Y is the scFv and when X is the scFv then Y is the peptide.

[0038] In one embodiment Y is independently selected from the 52SR4 scFv, and the GCN4 peptide, with the proviso that Y is the peptide X is the scFv and when Y is the scFv then X is the peptide.

[0039] In one embodiment the binding affinity between X and Y is 5 nM or stronger, for example the binding affinity of the heterodimeric tether is 900pM or stronger, such as 800, 700, 600, 500, 400 or 300pM.

[0040] In the present invention X or Y is a scFv 52SR4 (SEQ ID NO:3 or amino acids 1-243 of SEQ ID NO:3) specific to the peptide GCN4 (SEQ ID NO: 1 or amino acids 1 to 38 of SEQ ID NO: 1) and the other of X or Y is the peptide GCN4 (SEQ ID NO:1 or amino acids 1-38 of SEQ ID NO:1).

[0041] In one embodiment A and/or B is specific for an antigen selected from the group comprising: cell surface receptors such as T cell or B cell signalling receptors, co-stimulatory molecules, checkpoint inhibitors, natural killer cell receptors, Immunoglobulin receptors, immunoglobulin-like receptors, matrix metalloproteases and membrane type matrix metalloproteases tissue inhibitors of metalloproteases, TNFR family receptors, B7 family receptors, adhesion molecules, integrins, cytokine/chemokine receptors, GPCRs, growth factor receptors, kinase receptors, tissue-specific antigens, cancer antigens (tumour associated antigens & peptides), pathogen recognition receptors, complement receptors, hormone receptors, scavenger receptors, or soluble molecules such as cytokines, chemokines, leukotrienes, growth factors, hormones or enzymes or ion channels, including post translationally modified version thereof, fragments thereof comprising at least one epitope.

[0042] In one embodiment there is provided a composition, for example a pharmaceutical composition comprising one or more bispecific complexes according to the present disclosure.

[0043] Furthermore, the present disclosure also relates to a method that is not part of the invention of detecting synergistic function in a heterodimerically-tethered bispecific protein complex of the present invention, said method comprising the steps of:

(i) testing for activity in a functional assay for part or all of a multiplex comprising at least one heterodimerically-tethered bispecific protein; and

(ii) analysing the readouts from the functional assay to identify synergistic biological function in the bispecific protein complex.

[0044] The method employs the novel bispecific protein complex format of the present invention having the following formula A-X:Y-B wherein:

A-X is a first fusion protein;

Y-B is a second fusion protein;

X:Y is a heterodimeric-tether;

: is a binding interaction between X and Y;

A is a first protein component of the bispecific protein complex selected from a Fab or Fab' fragment;

B is a second protein component of the bispecific protein complex selected from a Fab or Fab' fragment;

wherein X is independently selected from a scFv, or a peptide, with the proviso that when X is a peptide Y is a scFv and when X is scFv then Y is a peptide;

wherein X or Y is a scFv specific to the peptide GCN4 represented by SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1; wherein the scFv is designated 52SR4 and is represented by SEQ ID NO:3 or amino acids 1 to 243 of SEQ ID NO:3 and wherein X or Y is a peptide GCN4 (SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1).

[0045] In particular, the heterodimerically-tethered bispecific protein complex is prepared by mixing A-X and B-Y *in vitro*. Thus in one embodiment the method comprises an *in vitro* mixing step bringing A-X and B-Y into contact.

[0046] Thus generally the fusion proteins A-X and B-Y are not co-expressed in the same cell. This is advantageous because it allows, for example 100 fusion proteins to be expressed and optionally purified and the subsequent mixing of the 100 fusion proteins in the various permutations can provide 10,000 heterodimerically-tethered bispecific protein complexes, of which 5,000 are unique pairs.

[0047] In contrast certain prior art methods require co-expression of bispecifics and thus for 10,000 complexes, 10,000 transfections, expressions and purifications are required.

[0048] However, if desired the A-X and B-Y may be expressed in the same cell.

[0049] The binding partners X and Y have affinity for each other and act as biological equivalent of velcro® or a bar and magnet and hold the complex together. Advantageously, this means that the fusion proteins A-X and Y-B can be readily assembled into a bispecific protein complex simply by mixing the fusion proteins together. Thus the bispecific protein complex of the present disclosure has a modular structure which allows for two different proteins to be easily

assembled in order to produce large panels of permutations of bispecific protein complexes with different combinations of antigen specificities in, for example a grid-like fashion. This allows for the efficient and systematic screening of a large number of bispecific protein complexes in order to detect additive, synergistic or novel biological function.

[0050] Given X and Y are specific for each other this significantly reduces the ability to form homodimers. X and Y are collectively referred to herein as a binding pair or binding partners. In one embodiment X does not have high affinity for other Xs. In one embodiment Y does not have high affinity for other Ys. Advantageously, when X and Y do not form homodimers, this prevents the formation of undesired monospecific protein complexes, increases yield of the desired bispecific protein complexes, and removes the need for onerous purification steps to remove the monospecific protein complexes.

[0051] This allows rapid assembly of bispecific protein complexes with a yield and/or purity which cannot be obtained efficiently by most prior art methods, in particular prior art methods generally require extensive purification steps. The yield of bispecific complex is typically 75% or higher in the present invention.

[0052] Further advantageously, the bispecific protein complexes allow for the screening of complexes wherein the constituent proteins (including antigens bound by the constituent proteins) do not have a known relationship or are in different potentially unrelated pathways, such as, two proteins which function in two distinct pathways and, for example which the skilled person would not normally expect to come into contact with each other can be tested in a bispecific protein complex to identify additive, synergistic and/or novel function.

[0053] Furthermore multiple binding regions (such as variable regions) to a given antigen or epitope can be investigated in parallel to identify nuances in biological function. This allows combinations of variable region sequences directed to a given pair of antigens to be investigated and optimised.

[0054] The present method allows the science to show the results and does not rely on pre-conceived ideas and technical prejudice about the biological function. This approach is potentially very powerful.

[0055] Advantageously the X and Y components allow a multiplex comprising bispecific protein complexes made up of different permutations of fusion proteins to be assembled rapidly and easily.

[0056] The proteins A and B are selected from Fab and Fab' fragments. When the Fab or Fab' fragments are held together as a complex via X and Y, this forms a bispecific antibody complex.

Description of Drawings

[0057]

Figure 1 is a schematic diagram showing the main cell signalling pathways involved in the activation of T cells.

Figure 2 is a schematic diagram showing the structure and assembly of a bispecific protein complex of the present disclosure.

Figure 3 is a table showing an example 4x4 grid for functional screening using the bispecific antibody of the present invention. Using this grid, 16 different bispecific protein complexes can be assembled and efficiently screened for synergistic function.

Figure 4 is a schematic diagram showing a representative bispecific antibody complex of the present disclosure. The diagram depicts how two different Fab fragments come together to form a bispecific antibody complex via the high affinity interaction between the binding partners attached to the Fab fragments.

Figure 5 is a graph showing the results of a flow cytometry experiment which demonstrates that two Fab fragments which are specific for two different target antigens retain their specificities and are able to co-engage their corresponding target antigens simultaneously when the two Fab fragments are combined to form a bispecific antibody complex of the present disclosure. The results further demonstrate that reversing the attachment of the binding partners attached to the Fab fragments does not affect the ability of the Fab fragments to bind specifically to their respective target antigens. The no complex formation control shows no binding detected when both specificities are fused to peptide (Y:Y). Filled area = [anti-antigen 5 Fab-peptide 'GCN4'] : [anti-antigen 6 Fab-peptide 'GCN4'] : [biotinylated-antigen 6] : [FITC-STREP] . No complex formation control. Thin line = [anti-antigen 5 Fab-scFv '52SR4'] : [anti-antigen 6 Fab-peptide 'GCN4'] : [biotinylated-antigen 6] : [FITC-STREP] Thick line = [anti-antigen 5 Fab-peptide 'GCN4'] : [anti-antigen 6 Fab-scFv '52SR4'] : [biotinylated-antigen 6] : [FITC-STREP]

Figure 6 is a graph and a table showing a BIAcore trace, which demonstrates the high affinity of the binding partners for each other. Fab A-scFv '52SR4' binding is detected to peptide 'GCN4' on the chip.

	Figure 7	is a bar chart of the relative potency of inhibition of phosphorylated Akt for bispecific and bivalent combinations of antibodies with specificity for antigen 3, antigen 1, antigen 4 and antigen 2.
5	Figure 8	is a bar chart of the relative potency of inhibition of phosphorylated PLCg2 for bispecific and bivalent combinations of antibodies with specificity for antigen 3, antigen 1, antigen 4 and antigen 2.
	Figure 9	is a bar chart of the relative potency of inhibition of CD86 expression for bispecific and bivalent combinations of antibodies with specificity for antigen 3, antigen 1, antigen 4 and antigen 2.
10	Figure 10	is a bar chart of the relative potency of inhibition of phosphorylated Akt for bispecific, bivalent or mixtures of antibodies with specificity for antigen 1 and antigen 2 as well as single Fab' controls.
	Figure 11	is a bar chart of the relative potency of inhibition of phosphorylated Akt for bispecific, bivalent or mixtures of antibodies with specificity for antigen 3 and antigen 2.
15	Figure 12	is a bar chart of the relative potency of inhibition of phosphorylated PLCg2 for bispecific, bivalent or mixtures of antibodies with specificity for antigen 3 and antigen 2.
	Figure 13	is a graph showing the titration of the effect of the bispecific combination of anti-antigen 3 and anti-antigen 2 on total I κ B levels in anti-IgM stimulated B cells.
20	Figure 14	is a graph showing the titration of the effect of the bispecific combination of antigen 3 and antigen 2 on CD86 expression on anti-IgM stimulated B cells.
	Figure 15	is a bar chart of the relative potency of inhibition of phosphorylated Akt for bispecific, bivalent or mixtures of antibodies with specificity for antigen 4 and antigen 2.
	Figure 16	is a bar chart of the relative potency of inhibition of phosphorylated PLCg2 for bispecific, bivalent or mixtures of antibodies with specificity for antigen 4 and antigen 2.
25	Figure 17	is a graph showing the titration of the effect of the bispecific combination of antigen 4 and antigen 2 on CD86 expression on anti-IgM stimulated B cells.
	Figure 18	is a graph showing the overlaid size exclusion A280 signal traces from experiment 1 of Example 11. The traces shown are: Fab-X (VR4247) control, Fab-Y (VR4248) control and 1 to 1 molar ratio mixture at 500 μ g/ml of Fab-X (VR4247) and Fab-Y (VR4248) complex. Peaks were detected at an absorbance of 280nm.
30	Figure 19	is a graph showing the overlaid size exclusion A214 signal traces from experiment 2 of Example 11. The traces shown are: Fab-X (VR4130) control, Fab-Y (VR4131) control and 1 to 1 molar ratio mixture at 500 μ g/ml of Fab-X (VR4130) and Fab-Y (VR4131) complex. Peaks were detected at an absorbance of 214nm.
35	Figure 20	is a graph showing the overlaid size exclusion A214 signal traces from experiment 2 of Example 11. The traces shown are all Fab-X (VR4130)/Fab-Y (VR4131) 1 to 1 molar ratio mixtures at 500 μ g/ml, 50 μ g/ml and 5 μ g/ml as indicated. Peaks were detected at an absorbance of 214nm.
40	Figure 21	is a table showing the data for the antigen grid cross specificities. Values are percentage inhibition (negative value for activation) of phosphorylation of Syk & represent the mean of multiple V region combinations evaluated.
	Figure 22	is a table showing the data for the antigen grid cross specificities. Values are percentage inhibition (negative value for activation) of phosphorylation of PLCg2 & represent the mean of multiple V-region combinations evaluated.
45	Figure 23	is a table showing the data for the antigen grid cross specificities. Values are percentage inhibition (negative value for activation) of phosphorylation of AKT & represent the mean of multiple V region combinations evaluated.
	Figure 24	is a graph showing the percentage inhibition of the phosphorylation of Syk, PLCg2 & AKT for each V-region combination for antigen 2 in Fab-X combined with antigen 3 in Fab-Y
50	Figure 25	is a graph showing the percentage inhibition of the phosphorylation of Syk, PLCg2 & AKT of the phosphorylation of Syk, PLCg2 & AKT for each V-region combination for antigen 3 in Fab-X combined with antigen 2 in Fab-Y.
	Figure 26	is a graph showing the percentage inhibition of the phosphorylation of Syk, PLCg2 & AKT for each V region combination for antigen 2 in Fab-X combined with antigen 4 in Fab-Y.
55	Figure 27	is a graph showing the percentage inhibition of the phosphorylation of Syk, PLCg2 & AKT for each V region combination for antigen 4 in Fab-X combined with antigen 2 in Fab-Y.
	Figure 28	is a graph showing the data for the percentage inhibition of anti-IgM induced CD71 expression on B-cells, by antigen3Fab'-X and antigen2-Fab'-Y when combination as either purified Fab'

- or from transient supernatant. **Circles** -Purified Antigen2Fab-Y + Antigen3Fab-X IC₅₀ 0.3224 nM **Squares** -Transient Sup Antigen2-Y + Antigen3-X IC₅₀ 0.2640 nM **Triangle** -Mock transfected supernatant control
- Figure 29** is a graph showing the data for the percentage inhibition of anti-IgM induced phosphorylation of p38 in B-cells, by antigen3-Fab'-X and antigen 2-Fab'-Y when combination as either purified Fab' or from transient supernatant. **Circles** -Purified Antigen2Fab-Y + Antigen3Fab-X IC₅₀ 0.1413 nM **Squares** -Transient Sup Antigen 2-Y + Antigen 3-X IC₅₀ 0.1861nM **Triangle** -Mock transfected supernatant control
- Key for Figure 30 to 33** 1. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6066); 2. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6078); 3. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6079); 4. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6080); 5. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6082); 6. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6067); 7. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6068); 8. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6070); 9. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6071); 10. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6073); 11. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6075); 12. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6076); 13. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6077); 14. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6069); 15. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6072); 16. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6074); 17. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (VR6081); 18. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (TSUP-24117); 19. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (TSUP-24432); 20. Mock supe 1; 21. Mock supe 2; 22. Antigen2Fab-Y (VR4447) + Antigen3Fab-X (4126) purified
- Figure 30** is a graph showing the inhibition of phospho readouts by purified antigen2-specific Fab-Y (VR4447) combined with antigen3-specific Fab-X transients on IgM stimulated B-cells from UCB_Cone_172.
- Figure 31** is a graph showing the inhibition of phosphor readouts by purified antigen 2-specific Fab-Y (VR4447) combined with antigen 3-specific Fab-X transients on IgM stimulated B-cells from UCB_Cone_173
- Figure 32** Inhibition of phospho readouts by purified antigen 2-specific Fab-Y (VR4450) combined with antigen 3-specific Fab-X transients on IgM stimulated B-cells from UCB_Cone_172
- Figure 33** Inhibition of phospho readouts by purified antigen 2-specific Fab-Y (VR4450) combined with antigen 3-specific Fab-X transients on IgM stimulated B-cells from UCB_Cone_173
- Figure 34** shows data for the percentage inhibition of anti-IgM induced phosphorylated PLCy2 in B-cells, by antigen 3 and antigen 2 specific Fab-Kd-Fab or BYbe
- Figure 35** shows data for the percentage inhibition of anti-IgM induced phosphorylated P38 in B-cells, by antigen 3 and antigen 2 specific Fab-Kd-Fab or BYbe
- Figure 36** shows data for the percentage inhibition of anti-IgM induced phosphorylated Akt in B-cells, by antigen 3 and antigen 2 specific Fab-Kd-Fab or BYbe
- Figure 37** shows data for the percentage inhibition of anti-IgM induced CD71 expression on B-cells, by antigen 3 and antigen 2 specific Fab-Kd-Fab or BYbe
- Figure 38** shows data for the percentage inhibition of anti-IgM induced CD40 expression on B-cells, by antigen 3 and antigen 2 specific Fab-Kd-Fab or BYbe.
- Figure 39** shows data for the percentage inhibition of anti-IgM induced CD86 expression on B-cells, by antigen 3 and antigen 2 specific Fab-Kd-Fab or BYbe
- Figure 40** shows the inhibition of CD27 expression on B cells by VR4447/VR4126 BYbe and VR4447/VR4126/VR645 BYbe/Albumin
- Figure 41** shows the inhibition of CD71 expression on B cells by VR4447/VR4126 BYbe and VR4447/VR4126/VR645 BYbe/Albumin
- Figure 42** shows the inhibition of CD86 expression on B cells by VR4447/VR4126 BYbe and VR4447/VR4126/VR645 BYbe/Albumin
- Figure 43** shows the inhibition of CD27 expression on B cells by VR4447/VR4130 BYbe and VR4447/VR4130/VR645 BYbe/Albumin
- Figure 44** shows the inhibition of CD71 expression on B cells by VR4447/VR4130 BYbe and VR4447/VR4130/VR645 BYbe/Albumin
- Figure 45** shows the inhibition of CD86 expression on B cells by VR4447/VR4130 BYbe and VR4447/VR4130/VR645 BYbe/Albumin

Detailed Description

[0058] "Bispecific protein complex" as used herein refers to a molecule comprising two proteins (A and B referred to herein as bispecific components also referred to herein as the first protein component and second protein component, respectively of the bispecific) which are retained together by a heterodimeric-tether. In one embodiment one or both of the proteins have a binding domain, for example one or both of the proteins are antibodies or fragments thereof (in particular a Fab or Fab' fragment, such complexes are also referred to as Fab-Kd-Fab).

[0059] "Fusion proteins" as employed herein comprise a protein component A or B fused to a binding partner X or Y (as appropriate). In one embodiment the fusion protein is a translational protein expressed by recombinant techniques from a genetic construct, for example expressed in a host from a DNA construct. In the context of the present disclosure one of the key characteristics of a fusion protein is that it can be expressed as a "single protein/unit" from a cell (of course in the case of fusion proteins comprising a Fab/Fab' fragment there will be two chains but this will be considered a single protein for the purpose of the present specification with one chain, typically the heavy chain fused at its C-terminus to X or Y as appropriate, optionally via a linker as described herein below).

[0060] The function of the heterodimeric tether X:Y is to retain the proteins A and B in proximity to each other so that synergistic function of A and B can be effected or identified, for example employing the method described herein.

[0061] "heterodimeric-tether" as used herein refers to a tether comprising two different binding partners X and Y which form an interaction : (such as a binding) between each other which has an overall affinity that is sufficient to retain the two binding partners together. In one embodiment X and/or Y are unsuitable for forming homodimers.

[0062] Heterodimerically-tethered and heterodimeric-tether are used interchangeably herein.

[0063] In one embodiment "unsuitable for forming homodimers" as employed herein refers to formation of the heterodimers of X-Y are more preferable, for example more stable, such as thermodynamically stable, once formed than homodimers. In one embodiment the binding interaction between X and Y is monovalent.

[0064] In one embodiment the X-Y interaction is more favourable than the X-X or Y-Y interaction.

[0065] This reduces the formation of homodimers X-X or Y-Y when the fusion proteins A-X and B-Y are mixed. Typically greater than 75% heterodimer is formed following 1:1 molar ratio mixing.

[0066] If desired, a purification step (in particular a one-step purification), such as column chromatography may be employed, for example to purify the fusion protein units and/or bispecific protein complexes according to the present disclosure.

[0067] In one embodiment a purification step is provided after expression of the or each fusion protein, although typically aggregate levels are low. Thus in one embodiment prior to *in vitro* mixing, the fusion protein(s) is/are provided in substantially pure form. Substantially pure form as employed herein refers to wherein the fusion protein is 90, 91, 92, 93, 94, 95, 96, 97, 98, 99 or 100% monomer.

[0068] In one embodiment no purification of the fusion protein or proteins is performed.

[0069] In one embodiment each fusion protein unit is expressed in a different expression experiment/run.

[0070] In one embodiment no purification of the fusion protein or proteins is performed before mixing to generate a bispecific protein complex. In one embodiment no purification of the fusion protein or proteins is performed before and/or after mixing.

[0071] In one embodiment no purification is required after the bispecific protein complex formation.

[0072] In one embodiment after mixing, and generally without further purification, at least 50% of the composition is the desired bispecific protein complex, for example at least 60, 65, 70, 75, 80% of the composition is the required bispecific protein complex.

[0073] In one embodiment the ratio of fusion proteins employed in the *in vitro* mixing step of the present method is A-X to B-Y 0.8:1 to 3:1, such as 1.5:1 or 2:1.

[0074] In one embodiment the ratio of fusion proteins employed in the *in vitro* mixing step of the present method is B-Y to A-X 0.8:1 to 3:1, such as 1.5:1 or 2:1, in a particular a molar ratio.

[0075] In one embodiment the ratio of A-X to B-Y employed in the *in vitro* mixing step is 1:1, in particular a 1:1 molar ratio.

[0076] The present disclosure also extends to a method of preparing a bispecific complex according to the present disclosure comprising admixing a fusion protein A-X and B-Y, for example in a 1:1 molar ratio.

[0077] In one embodiment the mixing occurs *in vitro*.

[0078] In one embodiment mixing occurs in a cell, for example a host cell.

[0079] In one embodiment, the mixing occurs *in vivo*, i.e. the fusion proteins A-X and B-Y interact with each other within a subject's body to form the heterodimeric-tether and in consequence, the bispecific protein complex.

[0080] In one embodiment, X and Y are completely specific for each other and do not bind to any other peptides/proteins in a cell or within a subject's body. This can be achieved for example by ensuring that X and Y are not naturally present in the target cell or in the target subject's body. This can be achieved, for example by selecting X or Y to be from a species or entity which is different to the subject (e.g. a yeast protein) and ensuring the other variable is specific to it. Advantageously, this prevents the binding of the fusion proteins A-X and/or B-Y to an undesired target, thereby generating

unwanted off-target effects.

[0081] In one embodiment one (or at least one) of the binding partners is incapable of forming a homodimer, for example an amino acid sequence of the binding partner is mutated to eliminate or minimise the formation of homodimers.

[0082] In one embodiment both of the binding partners are incapable of forming a homodimer.

[0083] Incapable of forming homodimers or aggregates as employed herein, refers to a low or zero propensity to form homodimers or aggregate. Low as employed herein refers to 5% or less, such as 4, 3, 2, 1, 0.5% or less aggregate, for example after mixing or expression or purification.

[0084] Small amounts of aggregate in the fusion proteins or residual in the heterodimerically-tethered bispecific protein complex generally has minimal effect on the screening method of the present disclosure. Therefore, in one embodiment no purification of fusion protein(s) and/or bispecific protein complex(es) is/are employed in the method, in particular after the mixing step.

[0085] In one embodiment : is a binding interaction based on attractive forces, for example Van der Waals forces, such as hydrogen bonding and electrostatic interactions, in particular, based on antibody specificity for an antigen (such as a peptide).

[0086] In one embodiment : is a covalent bond formed from a specific chemical interaction, such as click chemistry. In one embodiment : is not a covalent bond. In one embodiment conjugation/coupling chemistry is not employed to prepare the bispecific protein complexes of the present disclosure.

[0087] "Form the complex" as employed herein refers to an interaction, including a binding interaction or a chemical reaction, which is sufficiently specific and strong when the fusion protein components A-X and B-Y are brought into contact under appropriate conditions that the complex is assembled and the fusion proteins are retained together.

[0088] "Retained together" as employed herein refers to the holding of the components (the fusion proteins) in the proximity of each other, such that after X:Y binding the complex can be handled as if it were one molecule, and in many instances behaves and acts like a single molecule. In one embodiment the retention renders the complex suitable for use in the method disclosed herein, i.e. suitable for use in at least one functional screen.

[0089] Specificity as employed herein refers to where, for example the partners in the interaction e.g. X:Y or A and antigen or B and antigen only recognise each other or have significantly higher affinity for each other in comparison to non-partners, for example at least 2, 3, 4, 5, 6, 7, 8, 9, 10 times higher affinity, than for example a background level of binding to an unrelated non partner protein.

[0090] Specificity in relation to X and Y as employed herein refers to where the binding partners X and Y in the interaction only recognise each other or have significantly higher affinity for each other in comparison to non-partners, for example at least 2, 3, 4, 5, 6, 7, 8, 9, 10 times higher affinity.

[0091] In one embodiment the binding interaction is reversible. In one embodiment the binding interaction is essentially irreversible.

[0092] Essentially irreversible as employed herein refers to a slow off rate (dissociation constant) of the antibody or binding fragment.

[0093] In one embodiment, the binding interaction between X and Y has a low dissociation constant.

[0094] Examples of a low dissociation constant include $1-9 \times 10^{-2} \text{s}^{-1}$ or less, for example $1-9 \times 10^{-3} \text{s}^{-1}$, $1-9 \times 10^{-4} \text{s}^{-1}$, $1-9 \times 10^{-5} \text{s}^{-1}$, $1-9 \times 10^{-6} \text{s}^{-1}$ or $1-9 \times 10^{-7} \text{s}^{-1}$. Particularly suitable dissociation constants include $2 \times 10^{-4} \text{s}^{-1}$ or less, for example $1 \times 10^{-5} \text{s}^{-1}$, $1 \times 10^{-6} \text{s}^{-1}$ or $1 \times 10^{-7} \text{s}^{-1}$.

[0095] Whilst not wishing to be bound by theory it is thought that the low dissociation constant (also referred to as off rate) allows the molecules to be sufficiently stable to render the bispecific protein complex useful, in particular in functional screening assays.

[0096] In one embodiment, the affinity of X and Y for each other is 5 nM or stronger, for example 4nM, 3nM, 2nM, 1nM or stronger.

[0097] In one embodiment, the affinity of X and Y for each other is 900 pM or stronger, such as 800, 700, 600, 500, 400, 300, 200, 100 or 50 pM or stronger.

[0098] In another embodiment, the affinity of X and Y for each other is 10 pM or stronger, for example 9, 8, 7, 6 or 5 pM.

[0099] Affinity is a value calculated from the on and off rate of the entity. The term "affinity" as used herein refers to the strength of the sum total of non-covalent interactions between a single binding site of a molecule (e.g. an antibody) and its binding partner (e.g. a peptide).

[0100] The affinity of a molecule for its binding partner can generally be represented by the dissociation constant (KD). Affinity can be measured by common methods known in the art, including those described herein, such as surface plasmon resonance methods, in particular BIAcore.

[0101] However, the ability to hold the complex together is not just about affinity. Whilst not wishing to be bound by theory, we hypothesise that in fact there are three significant components: the on-rate, off-rate and the affinity. The calculation for affinity is based on on-rate and off-rate. So if the on-rate is low and the off-rate is fast, then the affinity will be low and that will not be sufficient to hold the bispecific protein complex together. However, a slow on-rate could be compensated for by a slow off-rate giving an overall suitable affinity. In some embodiments a high on-rate may be

sufficient to hold the complex together.

[0102] If the binding partners (X and Y) employed in the complex have a slow on-rate then additional time may be required after mixing the components to allow the complex to form.

[0103] If the affinity between the binding partners is sufficiently high, it may be possible for the bispecific protein complex to perform its desired biological function even if the affinity of the proteins (A and B) of the bispecific protein complex only bind weakly to their targets.

[0104] Conversely, if the proteins (A and B) are able to bind strongly to their targets, it may be possible to achieve the same biological function even if the affinity of the binding partners (X and Y) for each other is lower. In other words, a 'trinity' relationship exists such that a higher affinity between the binding partners can compensate for a lower affinity for the targets and *vice versa*.

[0105] In one embodiment the affinity of protein A for its ligand or antigen is about 100nM or stronger such as about 50nM, 20nM, 10nM, 1nM, 500pM, 250pM, 200pM, 100pM or stronger, in particular a binding affinity of 50pM or stronger.

[0106] In one embodiment the affinity of protein B for its ligand or antigen is about 100nM or stronger such as about 50nM, 20nM, 10nM, 1nM, 500pM, 250pM, 200pM, 100pM or stronger, in particular a binding affinity of 50pM or stronger.

[0107] In one embodiment an interaction between a constant domain in a heavy chain, such as CH1 and a constant domain in a light chain, such as CKappa contribute to the formation and/or stability of a bispecific complex according to the present disclosure. Thus employing Fab or Fab' fragments in the bispecific complexes of the present disclosure is beneficial.

[0108] In one embodiment the bispecific complex of the present disclosure does not comprise a component with an effector function, for example the complex does not comprise a constant domain other than a CH1 and CKappa or CLambda, in particular does not comprise constant domains independently selected from the group comprising CH2, CH3, CH4 and combinations thereof. In one embodiment the bispecific complex of the present disclosure lacks an Fc region.

[0109] In one embodiment the method herein is employed to screen a naive phage library by preparing fusion proteins of the disclosure from the library.

[0110] The bispecific protein complexes of the present invention may be used in any suitable application, including functional screening. This novel format is particularly useful in multiplex functional screening to identify protein targets based on function, and optimal epitopes on those target proteins, which could be targeted by bispecific therapies.

[0111] Furthermore where proteins A and B are antibodies or binding fragments thereof the bispecific protein complexes may also be used for multiplex functional screening to identify optimal variable region pairs for use in bispecific antibody therapeutics.

[0112] "Multiplex" as employed herein is a population of entities for testing comprising:

at least two component fusion proteins (A-X and Y-B) combined to generate at least one heterodimerically-tethered bispecific protein complex and at least one relevant biological comparator in the same or a different format, or at least two heterodimerically-tethered bispecific protein complexes with optionally at least one relevant biological comparator in the same or a different format.

[0113] Clearly to be useful, the different format employed as the comparator must be suitable for testing in a functional *in vitro* assay employed in the disclosure. In one example the comparator in the multiplex is a monovalent mixture of A-X and B-X or a bivalent monospecific complex of A-X-Y-A.

[0114] In one embodiment the multiplex comprises 1 to hundreds of thousands of heterodimerically-tethered bispecific protein complexes, for example 2 to 500,000 of said complexes, such as 2 to 100,000 or 2 to 10,000, in particular generated from mixing in a grid 2 to 100s of first and second fusion proteins (A-X and B-Y). In one embodiment the multiplex comprises for example 2 to 1,000, such as 2 to 900, 2 to 800, 2 to 700, 2 to 600, 2 to 500, 2 to 400, 2 to 300, 2 to 200, 2 to 100, 2 to 90, 3 to 80, 4 to 70, 5 to 60, 6 to 50, 7 to 40, 8 to 30, 9 to 25, 10 to 20 or 15 bispecific protein complexes. See Figure 3 for an example of such a grid.

[0115] In one embodiment the number of heterodimerically-tethered bispecific proteins in this multiplex is n^2 where n is 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15 or more.

[0116] The multiplex may be in the form of an array, for example a microtitre plate, wherein each well of the microplate may contain a different bispecific protein complex. The bispecific protein complexes may be tethered to a solid substrate surface, for example attached to a bead, or they may be suspended in a liquid (e.g. a solution or media) form, for example within a well or within a droplet.

[0117] In one embodiment every 'A' in the multiplex is a different protein, preferably an antibody or binding fragment thereof that binds to a target antigen and every 'B' is a different protein preferably an antibody or binding fragment thereof that binds to a target antigen.

[0118] In one embodiment the multiplex is provided in a grid as discussed below, for example an 8x8, 16x16 or 16x20, which equates to 64, 256 or 320 samples respectively.

[0119] "Grid" as employed herein refers to a two dimensional plot or array where one variable, such as protein A (in A-X) is varied along one axis, such as the X-axis (horizontal axis) and another variable such as protein B (in B-Y) is varied along the other axis, such as the Y axis (vertical axis). This arrangement assists in systematically evaluating the various combinations (permutations) of the variables.

[0120] In one embodiment the multiplex is provided on 96 well plates and the samples analysed may be multiples thereof i.e. 96, 192, 384 etc.

[0121] Advantageously, a grid arrangement is particularly advantageous for efficiently screening the biological function of bispecific protein complexes according to the present disclosure.

[0122] Figure 3 shows an example of such a grid, whereby 4 first fusion proteins can be readily combined with 4 second fusion proteins to produce 16 bispecific protein complexes.

[0123] Other variations of a screening grid will be apparent to the skilled addressee, for example the first protein (A) in the first fusion protein (A-X) may be kept constant whilst the second protein (B) in the second fusion protein (B-X) is varied. This may be useful for quickly screening a large number of different second proteins for synergistic function with the preselected first protein.

[0124] In another embodiment, protein A is varied along one axis by changing the antibody variable regions of protein A such that each antibody variant is specific for the same antigen but has a different combination of variable regions. Protein B may either be kept constant or may also be varied in the same fashion or varied such that the antigen specificity changes (across or down the grid) for the B proteins.

[0125] Advantageously, such a screening grid may potentially allow for the detection of small differences in synergistic function when the bispecific protein complexes are specific for the same antigens but with different combinations of variable regions.

[0126] In one embodiment, a "common" first fusion protein (A-X) according to the present disclosure may be present within each well. A range of different second fusion proteins (BY) according to the present disclosure may then be dispensed into each well. Subsequently, the specific binding interaction of the two binding partners (X and Y) physically brings the two fusion proteins together to form the bispecific protein complexes. This results in a multiplex comprising bispecific protein complexes which all bind to a common first target antigen (bound by A) but are also capable of binding to a second target antigen (bound by B) which may be different for each bispecific protein complex.

[0127] In one embodiment the B-Y fusion proteins comprise different variable regions to the same target antigen to allow optimisation of the variable regions and/or epitopes of the given target antigen bound by B when combined with the variable regions in A-X.

[0128] "Common" first fusion protein as employed herein refers to fusions proteins wherein the A or B component thereof, bind the same proteins or epitope, in particular where the A or B component have complete identity in the common fusion protein i.e. the common first fusion protein always comprises the same variable region sequence.

[0129] The skilled person is also aware of different variations of the above, such that the desired specificities of the bispecific protein complexes at each position in the multiplex can be readily controlled. This allows for the efficient screening of different combinations of bispecific protein complexes when such multiplexes are used in functional assays. In one embodiment factorial design is employed to define the variables employed in the grid.

[0130] In one embodiment the method of the present disclosure is conducive to high-throughput analysis.

[0131] In one embodiment, multiple bispecific protein complexes are tested in parallel or essentially simultaneously.

[0132] Simultaneously as employed herein refers to the where the samples/molecules/complexes are analysed in the same analysis, for example in the same "run". This may be advantageous as generally the reagents employed for a given sample run will be the same batch, concentration, cell source etc and therefore have the same properties. Furthermore the environmental conditions under which the analysis is performed, such as temperature and humidity are likely to be similar.

[0133] In one embodiment simultaneously refers to concomitant analysis where the signal output is analysed by an instrument at essentially the same time. This signal may require deconvolution to interpret the results obtained.

[0134] Advantageously, testing multiple bispecific protein complexes allows for more efficient screening of a large number of bispecific protein complexes and the identification of new and interesting relationships.

[0135] In one embodiment, the multiple bispecific protein complexes are tested by using a multiplex as defined above and subjecting the same to one or more functional assays. Accordingly the present invention provides a method for detecting synergistic biological function in a heterodimerically-tethered bispecific protein complex of formula A-X:Y-B

wherein X:Y is a heterodimeric-tether

: is a binding interaction between X and Y,

A and B are protein components of the bispecific in the form of fusion proteins with X and Y respectively, said method comprising the steps of:

(i) testing for activity in a functional assay for part or all of a multiplex comprising at least one heterodimerically-

tethered bispecific protein complex; and

(ii) analysing the readout(s) from the functional assay to identify or detect synergistic biological function in the heterodimerically-tethered bispecific protein complex; and

wherein Y is an antigen and X is an antibody or binding fragment thereof specific to Y or X is an antigen and Y is an antibody or binding fragment thereof specific to X.

[0136] The term "biological function" as used herein refers to an activity that is natural to or the purpose of, the biological entity being tested, for example a natural activity of a cell, protein or similar. Ideally the presence of the biological function can be tested using an *in vitro* functional assay, including assays employing mammalian cells, such as living cells, such as B or T cells, or tissue *ex vivo*. Natural function as employed herein also includes aberrant function, such as functions associated with diseases, such as cancers.

[0137] A relevant "biological comparator" as employed herein refers to a suitable entity for assessing activity, in the same assay as that employed for the bispecific protein complex, to establish if there is any change or novel activity or function. Suitable comparators for AX:Y-B may include a purified protein (including recombinant proteins) in a natural form or presented in the same format as the bispecific e.g., where A and B are the same entity, such as A-X:Y-A or B-X:Y-B i.e. a bivalent monospecific complex. Alternatively the fusion protein A-X or B-Y in an uncomplexed form may be employed as a comparator alone or as an uncomplexed mixture such as A-X and B-X together or A-Y and B-Y together.

[0138] Alternatively, multiple comparators of different formats (in particular as described herein) may be employed. The person skilled in the art is able to identify and include a suitable control/comparator based on common general knowledge or information that is found in the literature.

[0139] The term "synergistic function" or "synergistic biological function" as used herein refers to a biological activity or level of biological activity or an effect on a biological function or activity that:

- is not observed with individual fusion protein components until a bispecific is employed (and may include activity observed with a combination of antibodies to the said antigens, which are not in an bispecific format, but in particular refers to activity only observed when the two binding domains are linked in a bispecific format) or
- higher or lower activity in comparison to the activity observed when the first and second proteins of a bispecific protein complex of the present disclosure are employed individually, for example activity which is only observed in a bispecific form.

[0140] Therefore, "synergistic" includes novel biological function or novel activity. Synergistic function as employed herein does not generally include simple targeting i.e. based only on binding but will generally involve some inhibition, activation, signalling or similar after binding.

[0141] Novel biological function or novel activity as employed herein refers to a biological function or activity which is not apparent or is absent until the two or more synergistic entities (protein A and protein B) are brought together (as a bispecific or otherwise) or a previously unidentified function.

[0142] Higher as employed herein refers to an increase in activity including an increase from zero e.g., some activity in the bispecific where the individual uncomplexed bispecific component or components has/have no activity in the relevant functional assay, also referred to herein as new activity or novel biological function. Higher as employed herein also includes a greater than additive function in the bispecific in a relevant functional assay in comparison to the individual uncomplexed bispecific components (tested alone or in combination with being linked), for example 10, 20, 30, 40, 50, 60, 70, 80, 90, 100, 200, 300% or more increase in a relevant activity.

[0143] In one embodiment the uncomplexed proteins together have the same activity as the bispecific and this activity or function was previously unknown. This is also a novel synergistic function in the context of the present specification.

[0144] In one embodiment the synergistic function is a higher function.

[0145] In one embodiment the synergistic function is a lower function.

[0146] Lower function as employed herein refers to where the bispecific in the relevant functional assay has less or no activity in comparison to the individual uncomplexed bispecific component (s) which has/have activity in the relevant functional assay, also referred to herein as new activity or novel biological function (such as a natural protein i.e. a recombinant isolated protein which is not in a fusion protein nor part of any other complex other than one in which occurs *in vivo*-including an active domain or fragment of said protein) analysed as an individual protein or analysed as a mixture of proteins under the same conditions, for example 10, 20, 30, 40, 50, 60, 70, 80, 90, 100, 200, 300% or more decrease in a relevant activity. Greater than 100% decrease in activity refers to a gain in positive activity in a different direction, for example where an entity is an agonist decrease in activity over 100% may render the entity an antagonist and *vice versa*.

[0147] In one embodiment the activity of the bispecific complex is lower than the sum of the known function of protein A and protein B.

[0148] In some embodiments the bispecific protein complexes of the present disclosure have simply additive biological function. Additive biological function as employed herein refers to function, which is the same as the sum of each of the components A and B individually, when tested under the same conditions. An additive function may be a novel function if the activity or function was previously unknown or unidentified.

[0149] Screening is performed using any suitable assay known in the art, depending on the desired function to be identified.

[0150] In one embodiment, the functional assay employed in a method of the present disclosure is an *in vitro* or *ex vivo* assay.

[0151] A "functional assay," as used herein, is an assay that can be used to determine one or more desired properties or activities of the bispecific protein complexes, antibody complexes or the mixture of antibodies subject to the assay conditions. Suitable functional assays may be binding assays, apoptosis assays, antibody-dependent cellular cytotoxicity (ADCC) assays, complement-dependent cytotoxicity (CDC) assays, inhibition of cell growth or proliferation (cytostatic effect) assays, cell-killing (cytotoxic effect) assays, cell-signalling assays, cytokine production assays, antibody production and isotype switching, cellular differentiation assays, colony forming assays, chemotaxis assays, cell adhesion assays, cell migration assays, cell cycle assays, metabolic assays (whole cell and organelle function), assays for measuring inhibition of binding of pathogen to target cell, assays to measure the secretion of vascular endothelial growth factor (VEGF) or other secreted molecules, assays for bacteriostasis, bactericidal activity, neutralization of viruses, assays to measure the attraction of components of the immune system to the site where antibodies are bound, including *in situ* hybridization methods, labeling methods, and the like.

[0152] In one embodiment *in vivo* assays, such as animal models, including mouse tumor models, models of autoimmune disease, virus-infected or bacteria-infected rodent or primate models, and the like, may be employed.

[0153] The skilled person is well able to select a suitable functional assay based on the target/proteins being investigated. However, the complexes may be subject to a panel of "standard" assays without preselecting assays thought to be relevant in an attempt identify new functionality.

[0154] In the context of bispecific antibody complexes, the efficacy of bispecific antibody complexes according to the present disclosure can be compared to individual antibodies or mixtures of antibodies (or fragments) in such models by methods generally known to one of ordinary skill in the art.

[0155] For example, the bispecific antibody complexes may be tested for the ability to inhibit proliferation, affect viability or metabolic activity of cells (for example with a stain such as allamar blue or by monitoring luminescence due to luciferase expressed by the cells), or cause apoptosis of cancer cells, which are biological functions that include properties other than binding to an antigen.

[0156] By choosing functional assays closely related to a particular disease of interest, the methods of the disclosure make it possible to identify potentially therapeutic antibodies that bind to known or unknown target molecules. It is thus possible to identify new target molecules and/or to directly identify potentially therapeutic antibodies using the methods of the disclosure. Advantageously, the present method is not limited to any particular assay(s) and provides the user with complete flexibility to select the most appropriate functional assay depending on the requirements.

[0157] When screening the bispecific antibody complexes for desired biological function, various strategies may be employed. For example, medium containing the antibodies can be directly screened for the biological activity. Alternatively, the antibodies can be bound to beads coated or to microtiter plates prior to screening for biological activity. Alternatively a fusion protein maybe purified via a His tag in a nickel capture purification step. Such strategies may increase local concentrations of the antibodies leading to clearer results from the functional assays.

[0158] The functional assays may be repeated a number of times as necessary with or without different samples of a particular bispecific antibody complex to enhance the reliability of the results. Various statistical tests known to the skilled person can be employed to identify statistically significant results and thus identify bispecific antibody complexes with biological functions.

[0159] When establishing a functional assay for screening the skilled person can set a suitable threshold over which an identified activity is deemed a 'hit'. Where more than one functional assay is used the threshold for each assay may be set at a suitable level to establish a manageable hit rate. In one example the hit rate may be 3-5%. In one example the criteria set when searching for pairs of antigens that inhibit B cell function may be at least 30% inhibition of at least two phospho-readouts in a B cell activation assay.

[0160] In the bispecific protein complexes of the present invention the following protein and peptide components may be used.

[0161] In the present invention, X is independently selected from a scFv, or a peptide, with the proviso that when X is a peptide Y is a scFv and when X is scFv then Y is a peptide; wherein X or Y is a scFv specific to the peptide GCN4 represented by SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1; wherein the scFv is designated 52SR4 and is represented by SEQ ID NO:3 or amino acids 1 to 243 of SEQ ID NO:3 and wherein X or Y is a peptide GCN4 (SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1).

[0162] Other peptides referred to include the group comprising Fos/Jun (human and murine Fos have a Uniprot number

P01100 and P01101 respectively and human and murine jun have a Uniprot number 05412 and 05627 respectively), HA-tag which correspond to amino acids 98 to 106 of human influenza hemagglutinin, polyhistidine (His), c-myc and FLAG. Other peptides are also contemplated as suitable for use in the present disclosure and particularly suitable peptides are affinity tags for protein purification because such peptides have a tendency to bind with high affinity to their respective binding partners.

[0163] In one embodiment the peptide is not E5B9.

[0164] The term "peptide" as used herein refers to a short polymer of amino acids linked by peptide bonds, wherein the peptide contains in the range of 2 to 100 amino acids, for example 5 to 99, such as 6 to 98, 7 to 97, 8 to 96 or 5 to 25. In one embodiment a peptide employed in the present disclosure is an amino acid sequence of 50 amino acid residues or less, for example 40, 30, 20, 10 or less. The peptides used in the present disclosure are of a sufficient length to be fit for purpose, for example if the peptide is a linker, it needs to be suitably long to allow the fragment which it links to perform its biological function; alternatively if the peptide is a binding partner, it must be capable of binding specifically to another entity such as an antibody.

[0165] In one embodiment, the other binding partner of the binding pair (the alternative first or second binding partner) is a protein.

[0166] Protein as employed herein refers to an amino acid sequence of 100 amino acids or more. In one embodiment a "protein" as employed herein refers to an amino acid sequence with a secondary or tertiary structure.

[0167] Polypeptide and protein are employed interchangeably herein. However, polypeptide will generally be a protein with a simple structure, for example little secondary and/or tertiary structure.

[0168] In one embodiment the distinction between a peptide and a protein is based on the presence or absence of secondary structure and/or tertiary structure, where a peptide has no secondary structure and amino acids with secondary structure and/or tertiary structure are considered a protein.

[0169] In one embodiment, the protein is an antibody or an antibody fragment.

[0170] The term "antibody" as used herein refers to an immunoglobulin molecule capable of specific binding to a target antigen, such as a carbohydrate, polynucleotide, lipid, polypeptide, peptide etc., via at least one antigen recognition site (also referred to as a binding site herein), located in the variable region of the immunoglobulin molecule.

[0171] As used herein "antibody molecule" includes antibodies and binding fragments thereof. "Antibody fragments" as employed herein refer to antibody binding fragments including but not limited to Fab, modified Fab, Fab', modified Fab', F(ab')₂, Fv, single domain antibodies, scFv, bi, tri or tetra-valent antibodies, Bis-scFv, diabodies, triabodies, tetra-bodies and epitope-binding fragments of any of the above (see for example Holliger and Hudson, 2005, Nature Biotech. 23(9): 1126-1136; Adair and Lawson, 2005, Drug Design Reviews - Online 2(3), 209-217). The methods for creating and manufacturing these antibody fragments are well known in the art (see for example Verma et al., 1998, Journal of Immunological Methods, 216:165-181). Other antibody fragments for use in the present disclosure include the Fab and Fab' fragments described in International patent applications WO05/003169, WO05/003170 and WO05/003171. Multi-valent antibodies may comprise multiple specificities e.g. bispecific or may be monospecific (see for example WO92/22853, WO05/113605, WO2009/040562 and WO2010/035012).

[0172] A "binding fragment" as employed herein refers to a fragment capable of binding a target peptide or antigen with sufficient affinity to characterise the fragment as specific for the peptide or antigen.

[0173] The term "Fab fragment" as used herein refers to an antibody fragment comprising a light chain fragment comprising a VL (variable light) domain and a constant domain of a light chain (CL), and a VH (variable heavy) domain and a first constant domain (CH1) of a heavy chain. In one example the heavy chain sequences of the Fab fragment "terminates" at the interchain cysteine of CH1. In one embodiment the Fab fragment employed in a fusion protein of the present disclosure, such as A-X and/or B-Y is monovalent.

[0174] A Fab' fragment as employed herein refers to a Fab fragment further comprising all or part of a hinge region. In one embodiment the Fab' fragment employed in a fusion protein of the present disclosure, such as A-X and/or B-Y is monovalent.

[0175] The term "single-chain Fv" or abbreviated as "scFv", as used herein refers to an antibody fragment that comprises VH and VL antibody domains linked (for example by a peptide linker) to form a single polypeptide chain. The constant regions of the heavy and light chain are omitted in this format. Single-chain Fv as employed herein includes disulfide stabilised versions thereof wherein in addition to the peptide linker a disulfide bond is present between the variable regions.

[0176] Disulfide stabilised scFv may eliminate the propensity of some variable regions to dynamically breath, which relates to variable regions separating and coming together again.

[0177] The term "single domain antibody" as used herein refers to an antibody fragment consisting of a single monomeric variable antibody domain. Examples of single domain antibodies include VH or VL or VHH.

[0178] In one embodiment the antibody binding fragment and/or the bispecific antibody complex does not comprise an Fc region. "Does not comprise an Fc region" as employed herein refers to the lower constant domains, such as CH2, CH3 and CH4 which are absent. However, constant domains such as CH1, CKappa/CLambda may be present.

[0179] In one embodiment, the antibody heavy chain comprises a CH1 domain and the antibody light chain comprises

a CL domain, either kappa or lambda.

[0180] In one embodiment, the antibody heavy chain comprises a CH1 domain, a CH2 domain and a CH₃ domain and the antibody light chain comprises a CL domain, either kappa or lambda.

[0181] The first protein, A, and second protein, B, of the bispecific protein complex is an Fab or Fab' fragment. Such a bispecific protein complex may be referred to as a bispecific antibody complex.

[0182] "Bispecific antibody complex" as employed herein refers to a bispecific protein complex comprising at least two antibody binding sites wherein the component antibodies, fragments or both are complexed together by a heterodimeric-tether.

[0183] A bispecific antibody complex usually refers to a molecule comprising at least two antigen binding sites, wherein the binding sites have non-identical specificity.

[0184] In one embodiment, the two proteins (Fab or Fab' fragments) target the same antigen, for example binding to two different epitopes on the same target antigen, also referred to herein as a biparatopic bispecific.

[0185] In another embodiment, the two proteins (Fab or Fab' fragments) may have different antigen specificities, for example binding to two different target antigens.

[0186] In yet another embodiment, the two proteins are identical, i.e. binding to the same epitope on the same target antigen and the complex is therefore monospecific.

[0187] In one embodiment each antibody or fragment employed in the bispecific antibody complex of the disclosure comprises one binding site i.e. each binding site is monovalent for each target antigen.

[0188] The Fab or Fab' fragments employed in the fusion proteins (A-X or B-Y) may be monospecific, or monovalent.

[0189] Advantageously, the use of two bispecific antibody or antibody fragments allows the bispecific antibody complex of the present disclosure to potentially be specific for up to 4 different antigens (i.e. the complex may be tetraspecific). This allows avidity type effects to be investigated.

[0190] The antibody fragment employed in the first fusion protein (A-X) is a Fab or Fab'.

[0191] The fragment employed in the second fusion protein (B-Y) is a Fab, or Fab'.

[0192] "Monospecific" as employed herein refers to the ability to bind only one target antigen.

[0193] "Monovalent" as employed herein refers to the antibody or antibody fragment having a single binding site and therefore only binding the target antigen only once.

[0194] The Fab or Fab' fragment employed in the first fusion protein (A-X) is monovalent and the Fab or Fab' fragment employed in the second fusion protein (B-Y) is multivalent.

[0195] In one embodiment A-X or B-Y is not a fusion protein comprising two scFvs one specific to the antigen CD33 and one specific to the antigen CD3 or alternatively a bispecific complex format specific to these two antigens.

[0196] In one embodiment the A-X or B-Y is not a fusion protein comprising a scFv (or alternatively another antibody format) specific to CD3 linked to a peptide E5B9.

[0197] A "binding domain or site" as employed herein is the part of the antibody that contacts the antigen/epitope and participates in a binding interaction therewith. In one embodiment the binding domain contains at least one variable domain or a derivative thereof, for example a pair of variable domains or derivatives thereof, such as a cognate pair of variable domains or a derivative thereof.

[0198] In one embodiment the binding domain comprises 3 CDRs, in particular where the binding domain is a domain antibody such as a VH, VL or VHH. In one embodiment the binding domain comprises two variable domains and 6 CDRs and a framework and together these elements contribute to the specificity of the binding interaction of the antibody or binding fragment with the antigen/epitope.

[0199] A "cognate pair" as employed herein refers to a heavy and light chain pair isolated from a host as a pre-formed couple. This definition does not include variable domains isolated from a library, wherein the original pairings from a host is not retained. Cognate pairs may be advantageous because they are often affinity matured in the host and therefore may have high affinity for the antigen to which they are specific.

[0200] A "derivative of a naturally occurring domain" as employed herein is intended to refer to where one, two, three, four or five amino acids in a naturally occurring sequence have been replaced or deleted, for example to optimize the properties of the domain such as by eliminating undesirable properties but wherein the characterizing feature(s) of the domain is/are retained. Examples of modifications are those to remove glycosylation sites, GPI anchors, or solvent exposed lysines. These modifications can be achieved by replacing the relevant amino acid residues with a conservative amino acid substitution.

[0201] In one embodiment, the bispecific antibody complexes of the present disclosure or antibody/fragment components thereof are processed to provide improved affinity for a target antigen or antigens. Such variants can be obtained by a number of affinity maturation protocols including mutating the CDRs (Yang et al., J. Mol. Biol., 254, 392-403, 1995), chain shuffling (Marks et al., Bio/Technology, 10, 779-783, 1992), use of mutator strains of E. coli (Low et al., J. Mol. Biol., 250, 359-368, 1996), DNA shuffling (Patten et al., Curr. Opin. Biotechnol., 8, 724-733, 1997), phage display (Thompson et al., J. Mol. Biol., 256, 77-88, 1996) and sexual PCR (Cramer et al., Nature, 391, 288-291, 1998). Vaughan et al. (supra) discusses these methods of affinity maturation.

[0202] In one embodiment, the first Fab or Fab' fragment (A) is specific to a first antigen and the second Fab or Fab' fragment (B) is specific to a second antigen, wherein the first and second antigens are different. Advantageously, the bispecific antibody complex may be specific for two different antigens. This presents the possibility of the antibody complex binding to two different antigens, each located on a different entity, thereby bringing the two entities into close physical proximity with each other.

[0203] Alternatively, the first Fab or Fab' fragment (A) may be specific for a first epitope and the second Fab or Fab' (B) may be specific for a second epitope, wherein the first and second epitopes are both on the same antigen. This can greatly enhance the avidity of the bispecific antibody complex for the antigen due to the multiple interactions between the antigen and bispecific antibody complex.

[0204] In one embodiment, the first antibody (A) or second antibody (B) or both the first and second antibody of a bispecific antibody complex of the present disclosure may be an IgG, optionally with an inactive or active Fc region.

[0205] Also referred to, but which does not form part of the invention, is the situation where, the first (A) or second (B) antibody fragment is selected from the group consisting of: a fragment antigen binding (Fab), a Fab', a single chain variable fragment (scFv) and a single domain antibody (sdAb), such as a VHH

[0206] In one embodiment, the first antibody/fragment (A), second antibody/fragment (B) or both the first and second antibody/fragment of the bispecific antibody complex of the present disclosure may be a Fab.

[0207] In one embodiment, the first antibody/fragment (A), second antibody/fragment (B) or both the first and second antibody/fragment of the bispecific antibody complex of the present disclosure may be a Fab'.

[0208] For convenience bispecific protein complexes of the present disclosure are referred to herein as A-X:Y-B. However, this nomenclature is not intended to limit how the fusion protein A-X and B-Y are designed because our experiments indicate that binding partners X and Y can be reversed i.e. A-Y and B-X without adversely impacting on the method. Thus A and B and X and Y are nominal labels referred to for assisting the explanation of the present technology.

[0209] "Attached" as employed herein refers to connected or joined directly or indirectly via a linker, such as a peptide linker examples of which are discussed below. Directly connected includes fused together (for example a peptide bond) or conjugated chemically.

[0210] "Binding partner" as employed herein refers to one component part of a binding pair.

[0211] In one embodiment, the affinity of the binding partners is high, 5 nM or stronger, such as 900, 800, 700, 600, 500, 400, 300 pM or stronger.

[0212] "Binding pair" as employed herein refers to two binding partners which specifically bind to each other. Examples of a binding pair include a peptide and an antibody or binding fragment specific thereto, or an enzyme and ligand, or an enzyme and an inhibitor of that enzyme.

[0213] The first binding partner (X) is the scFv 52SR4 of SEQ ID NO: 3 or amino acids 1 to 243 of SEQ ID NO:3 or is the peptide GCN4 of SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1.

[0214] When X is the scFv 52SR4 then Y is the peptide GCN4.

[0215] The second partner (Y) is the scFv 52SR4 of SEQ ID NO: 3 or amino acids 1 to 243 of SEQ ID NO:3 or is the peptide GCN4 of SEQ ID NO: 1 or amino acids 1 to 38 of SEQ ID NO:1. When Y is the scFv 52SR4 then X is the peptide GCN4.

[0216] In one embodiment, the first binding partner (X) is attached to the C-terminal of the heavy or light chain of the Fab or Fab' fragment of, for example, the first binding partner (X) is attached to the C-terminal of the heavy chain of the Fab or Fab' fragment (A).

[0217] In another embodiment, the second binding partner (Y) is attached to the C-terminal of the heavy or light chain of the second Fab or Fab' fragment of B, for example the second binding partner (Y) is attached to the C-terminal of the heavy chain of the Fab or Fab' fragment (B).

[0218] In one embodiment X is attached to the C-terminal of the heavy chain of the Fab or Fab' fragment (protein A) and Y is attached to the C-terminal of the heavy chain of the Fab or Fab' fragment (protein B).

[0219] In one embodiment X is attached via a linker (such as ASGGGG SEQ ID NO: 71 or ASGGGGSG SEQ ID NO: 72) or any other suitable linker known in the art or described herein below, to the C-terminal of the heavy chain of the antibody or fragment (protein A) and Y is attached via a linker (such as ASGGGG SEQ ID NO: 71 or ASGGGGSG SEQ ID NO: 72) to the C-terminal of the heavy chain of the antibody or fragment (protein B).

[0220] Examples of a suitable binding pair (X or Y) may include GCN4 (SEQ ID NO: 1 or lacking the HIS tag, amino acids 1-38 of SEQ ID NO: 1) or a variant thereof and 52SR4 (SEQ ID NO: 3 or lacking the HIS tag amino acids 1 to 243 of SEQ ID NO:3) or a variant thereof, which is a scFv specific for GCN4.

[0221] In a one embodiment of the invention the first binding partner (nominally X) is GCN4 (as shown in SEQ ID NO: 1) or a variant thereof without the His tag and the second binding partner (nominally Y) is a scFv specific for GCN4 (as shown in SEQ ID NO: 3 or a variant thereof represented by amino acids 1 to 243 of SEQ ID NO: 3).

[0222] In one embodiment of the invention, the first binding partner (nominally X) is a scFv specific for GCN4 (as shown in SEQ ID NO: 3) and the second binding partner (nominally Y) is GCN4 (as shown in SEQ ID NO: 1).

[0223] The present inventors have found that the single chain antibody 52SR4 and peptide GCN4, are a binding pair suitable for use in the bispecific protein complexes of the present disclosure.

[0224] "Active fragment" as employed herein refers to an amino acid fragment, which is less than the whole amino acid sequence for the entity and retains essentially the same biological activity or a relevant biological activity, for example greater than 50% activity such as 60%, 70%, 80%, 85%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100%.

[0225] "Identity", as used herein, indicates that at any particular position in the aligned sequences, the amino acid residue is identical between the sequences. "Similarity", as used herein, indicates that, at any particular position in the aligned sequences, the amino acid residue is of a similar type between the sequences. For example, leucine may be substituted for isoleucine or valine. Other amino acids which can often be substituted for one another include but are not limited to:

- phenylalanine, tyrosine and tryptophan (amino acids having aromatic side chains);
- lysine, arginine and histidine (amino acids having basic side chains);
- aspartate and glutamate (amino acids having acidic side chains);
- asparagine and glutamine (amino acids having amide side chains); and
- cysteine and methionine (amino acids having sulphur-containing side chains).

[0226] Degrees of identity and similarity can be readily calculated (Computational Molecular Biology, Lesk, A.M., ed., Oxford University Press, New York, 1988; Biocomputing. Informatics and Genome Projects, Smith, D.W., ed., Academic Press, New York, 1993; Computer Analysis of Sequence Data, Part 1, Griffin, A.M., and Griffin, H.G., eds., Humana Press, New Jersey, 1994; Sequence Analysis in Molecular Biology, von Heinje, G., Academic Press, 1987, Sequence Analysis Primer, Gribskov, M. and Devereux, J., eds., M Stockton Press, New York, 1991, the BLAST™ software available from NCBI (Altschul, S.F. et al., 1990, J. Mol. Biol. 215:403-410; Gish, W. & States, D.J. 1993, Nature Genet. 3:266-272. Madden, T.L. et al., 1996, Meth. Enzymol. 266:131-141; Altschul, S.F. et al., 1997, Nucleic Acids Res. 25:3389-3402; Zhang, J. & Madden, T.L. 1997, Genome Res. 7:649-656.).

[0227] In one embodiment, the first or second binding partner (X or Y) is a protein or peptide.

[0228] In one embodiment, the first and second fusion proteins comprise one or more peptide linkers. The linkers may be incorporated at various locations in the fusion proteins. For example, a linker may be introduced between a binding partner and the protein attached thereto.

[0229] In one embodiment, the linker is a peptide linker.

[0230] The term "peptide linker" as used herein refers to a peptide with an amino acid sequence. A range of suitable peptide linkers will be known to the person of skill in the art.

[0231] In one embodiment, the binding partners of the bispecific protein complexes are joined to their respective proteins via peptide linkers.

[0232] In one embodiment the fusion proteins are a translational fusion, that is a fusion protein expressed in a host cell comprising a genetic construct from which the fusion protein is expressed.

[0233] In one embodiment the fusion protein is prepared by fusing the heavy chain of A to X and/or the heavy chain of B to Y optionally via a peptide linker.

[0234] In one embodiment, the peptide linker is 50 amino acids in length or less, for example 20 amino acids or less.

[0235] Generally it will be more efficient to express the fusion protein recombinantly and therefore a direct peptide bond or a peptide linker that can be expressed by a host cell may be advantageous.

[0236] In one embodiment, the linker is selected from a sequence shown in sequence 5 to 72 or PPP.

Table 1. Hinge linker sequences

SEQ ID NO:	SEQUENCE
5	DKTHTCAA
6	DKTHTCPPCPA
7	DKTHTCPPCPATCPPCPA
8	DKTHTCPPCPATCPPCPATCPPCPA
9	DKTHTCPPCPAGKPTLYNSLVMSDAGTCY
10	DKTHTCPPCPAGKPTHVNVSVVMAEVDGTCY
11	DKTHTCCVECPPCPA
12	DKTHTCPRCPEPKSCDTPPPCPRCPA

(continued)

SEQ ID NO:	SEQUENCE
13	DKTHTCPSCPA

Table 2. Flexible linker sequences

SEQ ID NO:	SEQUENCE
14	SGGGGSE
15	DKTHTS
16	(S)GGGGS
17	(S)GGGGSGGGGS
18	(S)GGGGSGGGGS
19	(S)GGGGSGGGGS
20	(S)GGGGSGGGGS
21	AAAGSG-GASAS
22	AAAGSG-XGGGS-GASAS
23	AAAGSG-XGGGSXGGGS -GASAS
24	AAAGSG- XGGGSXGGGSXGGGS -GASAS
25	AAAGSG- XGGGSXGGGSXGGGSXGGGS-GASAS
26	AAAGSG-XS-GASAS
27	PGGNRGTTTTRRPATTTGSSPGPTQSHY
28	ATTTGSSPGPT
29	ATTTGS
30	GS
31	EPSGPISTINSPPSKESHKSP
32	GTVAAPSVFIFPPSD
33	GGGGIAPSMVGGGGS
34	GGGGKVEGAGGGGGS
35	GGGGSMKSHDGGGGS
36	GGGGLITIVGGGGS
37	GGGGVPSLPGGGGS
38	GGEKSIPGGGGS
39	RPLSYRPPFPFGFPSVRP
40	YPRSIYIRRRHPSPSLTT
41	TPSHLSHILPSFGLPTFN
42	RPVSPFTFPRLNSWLP A
43	SPAAHFPRSIPRPGPIRT
44	APGPSAPSHRSLPSRAFG
45	PRNSIHFLHPLLVAPLGA
46	MPSLSGVLQVRYLSPPDL

(continued)

SEQ ID NO:	SEQUENCE
47	SPQYPSPLTLTLPHPHPSL
48	NPSLNPPSYLHRAPSRIS
49	LPWRTSLLPSLPLRRRP
50	PPLFAKGPVGLLSRSFPP
51	VPPAPVVSLRSAHARPPY
52	LRPTPPRVRSYTCCPTP-
53	PNVAHVLPLLTVPWDNLR
54	CNPLLPLCARSPAVRTFP
(S) is optional in sequences 17 to 20.	

[0237] Examples of rigid linkers include the peptide sequences GAPAPAAPAPA (SEQ ID NO: 69), PPPP (SEQ ID NO: 70) and PPP.

[0238] Other linkers are shown in **Table 3**:

SEQ ID NO:	SEQUENCE
55	DLCLRDWGCLW
56	DICLPRWGCLW
57	MEDICLPRWGCLWGD
58	QRLMEDICLPRWGCLWEDDE
59	QGLIGDICLPRWGCLWGRSV
60	QGLIGDICLPRWGCLWGRSVK
61	EDICLPRWGCLWEDD
62	RLMEDICLPRWGCLWEDD
63	MEDICLPRWGCLWEDD
64	MEDICLPRWGCLWED
65	RLMEDICLARWGCLWEDD
66	EVRSFCTRWPAEKSKPLRG
67	RAPESFVCYWETICFERSEQ
68	EMCYFPGICWM

[0239] In one aspect of the disclosure that is not part of the invention there is provided a method of producing a bispecific protein complex of the present disclosure, comprising the steps of:

- (a) producing a first fusion protein (A-X), comprising a first protein (A), attached to a first binding partner (X) of a binding pair;
- (b) producing a second fusion protein (B-Y), comprising a second protein (B), attached to a second binding partner (Y) of a binding pair; and
- (c) mixing the first (A-X) and second fusion proteins (B-Y) prepared in step a) and b) together.

[0240] Typically the mixing of A-X and B-Y in step (c) is in a 1:1 molar ratio.

[0241] In one embodiment each fusion proteins employed in the complexes of the present disclosure are produced by expression in a host cell or host cells in an expression experiment.

[0242] In one aspect, there is provided a method of preparing a bispecific protein complex of the present disclosure,

comprising the steps of:

- (a) expressing a first fusion protein (A-X), comprising a first protein (A), attached to a first binding partner (X) of a binding pair;
- (b) expressing a second fusion protein (B-Y), comprising a second protein (B), attached to a second binding partner (Y) of a binding pair;

wherein fusion protein A-X and B-Y are expressed from the same host cell or distinct host cells. In said method the bispecific protein complex produced has the formula A-X:Y-B wherein:

A-X is a first fusion protein;
Y-B is a second fusion protein;
X:Y is a heterodimeric-tether;

wherein X is independently selected from a scFv, or a peptide, with the proviso that when X is a peptide Y is a scFv and when X is scFv then Y is a peptide;
wherein X or Y is a scFv specific to the peptide GCN4 (SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1); wherein the scFv is 52SR4 (SEQ ID NO:3 or amino acids 1 to 243 of SEQ ID NO:3) and wherein X or Y is a peptide GCN4 (SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO: 1).

[0243] Distinct host cells as employed herein refers to individual cells, including cells of the same type (even same clonal type).

[0244] In one embodiment the expression is transient expression. The use of transient expression is highly advantageous when combined with the ability to generate bispecific complexes without recourse to purification. This results in a rapid method to generate bispecific protein complexes as transient transfection is much simpler and less resource intensive than stable transfection.

[0245] In one embodiment the expression is stable expression i.e. wherein the DNA encoding the fusion protein in question is stably integrated into the host cell genome.

[0246] In one embodiment a polynucleotide encoding A-X and a polynucleotide encoding B-Y on the same or different polynucleotide sequences are transfected into a cell as part of a functional assay, wherein the proteins are expressed in the cell and/or released therefrom. In particular the polynucleotides are transiently transfected on the same of different plasmids.

[0247] The mixing of A-X and B-Y is generally effected in conditions where the X and Y can interact. In one embodiment, the fusion proteins are incubated in cell culture media under cell culturing conditions, for example the fusion proteins are incubated for 90 minutes in a 37°C/5%CO₂ environment.

[0248] In one embodiment the fusion proteins of the present disclosure are mixed in an aqueous environment, for example one fusion protein may be bound to a solid surface such as a bead or a plate and the other fusion protein can be introduced thereto in an aqueous solution/suspension. The solid phase allows excess components and reagents to be washed away readily. In one embodiment neither fusion is attached a solid phase and are simply mixed in a liquid/solution/medium. Thus in one embodiment A-X and B-Y are mixed as free proteins in an aqueous media.

[0249] Advantageously, the method of the present disclosure can be employed to prepare complexes formed between heterogenous pairs (i.e. between the first fusion protein [A-X] and second fusion protein [B-Y]) wherein interactions between homogenous pairs (i.e. between two first fusion proteins [A-X] or two second fusion proteins [B-Y]) are minimised. Thus the present method allows large numbers of bispecific protein complexes to be prepared, with minimal or no contamination with homodimeric complexes. An advantage of the constructs and method of the present disclosure is that the ratio of A-X to B-Y is controlled by the properties of the A-X and B-Y and in particular a molar ratio of 1:1 can be achieved. This element of control is a significant improvement over the certain prior art methods.

[0250] In one embodiment a method of the present disclosure comprises a further step of transferring a pair variable regions (in particular two pairs of variable regions) identified as having synergistic activity into an alternative bispecific format, optionally humanising said variable regions if necessary beforehand, which is an alternative therapeutic format and/or a format having an extended half-life suitable for testing in assays with a longer duration (for example which run a week or more).

[0251] Multivalent formats include those known in the art and those described herein, such as DVD-IgS, FabFvs for example as disclosed in WO2009/040562 and WO2010/035012, diabodies, triabodies, tetrabodies etc.

[0252] Other examples of bi and multispecific formats (including therapeutic formats) include a diabody, triabody, tetrabody, tandem scFv, tandem scFv-Fc, FabFv, Fab'Fv, FabdsFv, Fab-scFv, Fab'-scFv, diFab, diFab', scdiabody, scdiabody-Fc, ScFv-Fc-scFv, scdiabody-CH₃, IgG-scFv, scFv-IgG, V-IgG, IgG-V, DVD-Ig, and DuoBody.

[0253] Diabody as employed herein refers to two Fv pairs: VH/VL and a further VH/VL pair which have two inter-Fv linkers, such that the VH of a first Fv is linked to the VL of the second Fv and the VL of the first Fv is linked to the VH of

the second Fv.

[0254] Triabody as employed herein refers to a format similar to the diabody comprising three Fv pairs and three inter-Fv linkers.

[0255] Tetrabody as employed herein refers to a format similar to the diabody comprising four Fv pairs and four inter-Fv linkers.

[0256] Tandem scFv as employed herein refers to two scFvs (each comprising a linker in the usual manner) linked to each other via a single linker such that there is a single inter-Fv linker.

[0257] Tandem scFv-Fc as employed herein refers to two tandem scFvs, wherein each one is appended to the N-terminus of a CH2 domain, for example via a hinge, of constant region fragment -CH2CH3.

[0258] FabFv as employed herein refers to a Fab fragment with a variable region appended to the C-terminal of each of the following, the CH1 of the heavy chain and CL of the light chain. The format may be provided as a PEGylated version thereof.

[0259] Fab'Fv as employed herein is similar to FabFv, wherein the Fab portion is replaced by a Fab'. The format may be provided as a PEGylated version thereof.

[0260] FabdsFv as employed herein refers to a FabFv wherein an intra-Fv disulfide bond stabilises the appended C-terminal variable regions. The format may be provided as a PEGylated version thereof.

[0261] Fab-scFv as employed herein is a Fab molecule with a scFv appended on the C-terminal of the light or heavy chain.

[0262] Fab'-scFv as employed herein is a Fab' molecule with a scFv appended on the C-terminal of the light or heavy chain.

[0263] DiFab as employed herein refers to two Fab molecules linked via their C-terminus of the heavy chains.

[0264] DiFab' as employed herein refers to two Fab' molecules linked via one or more disulfide bonds in the hinge region thereof.

[0265] As employed herein scdiabody is a diabody comprising an intra-Fv linker, such that the molecule comprises three linkers and forms a normal scFv whose VH and VL terminals are each linked to a one of the variable regions of a further Fv pair.

[0266] Scdiabody-Fc as employed herein is two scdiabodies, wherein each one is appended to the N-terminus of a CH2 domain, for example via a hinge, of constant region fragment -CH2CH3.

[0267] ScFv-Fc-scFv as employed herein refers to four scFvs, wherein one of each is appended to the N-terminus and the C-terminus of both the heavy and light chain of a -CH2CH3 fragment.

[0268] Scdiabody-CH3 as employed herein refers to two scdiabody molecules each linked, for example via a hinge to a CH3 domain.

[0269] IgG-scFv as employed herein is a full length antibody with a scFv on the C-terminal of each of the heavy chains or each of the light chains.

[0270] scFv-IgG as employed herein is a full length antibody with a scFv on the N-terminal of each of the heavy chains or each of the light chains.

[0271] V-IgG as employed herein is a full length antibody with a variable domain on the N-terminal of each of the heavy chains or each of the light chains.

[0272] IgG-V as employed herein is a full length antibody with a variable domain on the C-terminal of each of the heavy chains or each of the light chains.

[0273] DVD-Ig (also known as dual V domain IgG) is a full length antibody with 4 additional variable domains, one on the N-terminus of each heavy and each light chain.

[0274] Duobody or 'Fab-arm exchange' as employed herein is a bispecific IgG antibody format where matched and complementary engineered amino acid changes in the constant domains (typically CH3) of two different monoclonal antibodies lead, upon mixing, to the formation of heterodimers. A heavy/light chain pair from the first antibody will, as a result of the residue engineering, prefer to associate with a heavy/light chain pair of a second antibody.

[0275] If present constant region domains of a bispecific antibody complex or antibody molecule of the present disclosure, if present, may be selected having regard to the proposed function of the complex or antibody molecule, and in particular the effector functions which may be required. For example, the constant region domains may be human IgA, IgD, IgE, IgG or IgM domains. In particular, human IgG constant region domains may be used, especially of the IgG1 and IgG3 isotypes when the antibody molecule is intended for therapeutic uses and antibody effector functions are required. Alternatively, IgG2 and IgG4 isotypes may be used when the antibody molecule is intended for therapeutic purposes and antibody effector functions are not required. It will be appreciated that sequence variants of these constant region domains may also be used. For example IgG4 molecules in which the serine at position 241 has been changed to proline as described in Angal et al., 1993, Molecular Immunology, 1993, 30:105-108 may be used. Accordingly, in the embodiment where the antibody is an IgG4 antibody, the antibody may include the mutation S241P.

[0276] It will also be understood by one skilled in the art that antibodies may undergo a variety of posttranslational modifications. The type and extent of these modifications often depends on the host cell line used to express the antibody

as well as the culture conditions. Such modifications may include variations in glycosylation, methionine oxidation, diketopiperazine formation, aspartate isomerization and asparagine deamidation. A frequent modification is the loss of a carboxy-terminal basic residue (such as lysine or arginine) due to the action of carboxypeptidases (as described in Harris, R.J. Journal of Chromatography 705:129-134, 1995). Accordingly, the C-terminal lysine of the antibody heavy chain may be absent.

[0277] The present disclosure also provides a composition comprising one or more bispecific protein complexes as described above, wherein the composition predominantly comprises heterodimeric bispecific complexes according to the present disclosure, for example with minimal or no contamination with homodimeric complexes.

[0278] In one embodiment, at least 60%, at least 65%, at least 70%, at least 75%, at least 80%, at least 90%, or at least 95% of the fusion proteins in the composition are in a bispecific protein complex form.

[0279] In one embodiment, at least 60% of the fusion proteins in the composition are in a bispecific protein complex form.

[0280] In one embodiment the complexes formed require no further purification steps and thus the compositions comprise unpurified bispecific complexes.

[0281] In one embodiment the complexes formed require one purification step, for example column chromatography.

[0282] In one embodiment the method further comprises at least one purification step, for example after expression of a fusion protein according to the present disclosure and before mixing the fusion proteins.

[0283] In one aspect the present disclosure relates to a fusion protein, a heterodimerically-tethered bispecific protein complex, a composition comprising a fusion protein or said bispecific protein complex, a multiple, array, library as defined herein.

[0284] In one embodiment, the bispecific protein complex is in solution or suspension.

[0285] In one embodiment, the bispecific protein complexes are fixed on a solid substrate surface.

[0286] In one embodiment, the multiplex is in the form of an array, for example in a microplate, such as a 96 or 384 well plate. Such arrays can be readily implemented in screening assays to identify bispecific protein complexes with desired functionality.

[0287] In another embodiment, the bispecific protein complexes are conjugated to beads.

[0288] A fusion protein as defined above is a component of the bispecific protein complex according to the present disclosure. In one aspect, the present disclosure relates to a fusion protein described herein.

[0289] In a further aspect, there is provided a library, comprising two or more fusion proteins as defined above.

[0290] The term "library" as used herein refers to two or more bispecific antibody complexes of the present disclosure or multiple fusion proteins of the present disclosure that can be combined to form at least two different bispecific antibody complexes according to the present disclosure. As described throughout the specification, the term "library" is used in its broadest sense and may also encompass sub-libraries.

[0291] Advantageously, the library may comprise a range of different fusion proteins which have either the first binding partner (X) or second binding partner (Y) of a particular binding pair attached thereto. In one embodiment part of the library comprises proteins/antibodies/fragments each connected to a binding partner X and the remainder of the library comprises the same proteins/antibodies/fragments each connected to a binding partner Y. This thus allows any two fusion proteins to be readily combined to form a bispecific protein complex of the present disclosure, as long as one fusion protein has the first binding partner of a binding pair attached and the other fusion protein has the second binding partner of the binding pair attached.

[0292] In one embodiment bispecific protein complexes of the present invention are suitable for therapeutic applications and may provide novel therapies for treating diseases. Thus in a further aspect, there is provided a bispecific protein complex as described above for use in therapy. The bispecific protein complex is suitable for treating a range of diseases, such as autoimmune disease and cancer.

[0293] Conversely, the bispecific protein complexes of the present disclosure can be engineered with one antibody or antibody fragment specific for T-lymphocytes, and another antibody or antibody fragment specific for a cancer-specific antigen. As a result, the bispecific antibody complexes of the present disclosure may advantageously possess a higher cytotoxic potential compared to ordinary monoclonal antibodies.

[0294] The bispecific protein complexes of the present disclosure are also particularly suited for inhibiting B cell function in order to control immune and autoimmune reactions in various autoimmune diseases.

[0295] Thus, the present disclosure extends to a method of treating a disease in a patient, comprising the administration of a bispecific protein complex of the present disclosure.

[0296] In one aspect, there is provided a pharmaceutical composition comprising one or more bispecific protein complexes of the present disclosure.

[0297] In one embodiment there is provided a fusion protein obtained or obtainable for a method of the present disclosure.

[0298] In one embodiment there is provided an bispecific antibody complex obtained or obtainable from a method of the present disclosure

[0299] In one embodiment there is provided an a bispecific or multispecific antibody molecule comprising variable

regions combinations identified by a method according to the present disclosure.

[0300] In one embodiment there is provided a composition, such as a pharmaceutical composition comprising a fusion protein, a bispecific antibody complex or a bispecific/multispecific antibody molecule obtained from a method of the present disclosure.

[0301] Various different components can be included in the composition, including pharmaceutically acceptable carriers, excipients and/or diluents. The composition may, optionally, comprise further molecules capable of altering the characteristics of the population of antibodies of the invention thereby, for example, reducing, stabilizing, delaying, modulating and/or activating the function of the antibodies. The composition may be in solid, or liquid form and may *inter alia*, be in the form of a powder, a tablet, a solution or an aerosol.

[0302] The present disclosure also provides a pharmaceutical or diagnostic composition comprising a bispecific protein complex of the present invention in combination with one or more of a pharmaceutically acceptable excipient, diluent or carrier. Accordingly, provided is the use of a bispecific protein complex of the invention for use in the treatment and for the manufacture of a medicament for the treatment of a pathological condition or disorder.

[0303] The pathological condition or disorder, may, for example be selected from the group consisting of infections (viral, bacterial, fungal and parasitic), endotoxic shock associated with infection, arthritis such as rheumatoid arthritis, asthma such as severe asthma, chronic obstructive pulmonary disease (COPD), pelvic inflammatory disease, Alzheimer's Disease, inflammatory bowel disease, Crohn's disease, ulcerative colitis, Peyronie's Disease, coeliac disease, gallbladder disease, Pilonidal disease, peritonitis, psoriasis, vasculitis, surgical adhesions, stroke, Type I Diabetes, Lyme disease, meningoencephalitis, autoimmune uveitis, immune mediated inflammatory disorders of the central and peripheral nervous system such as multiple sclerosis, lupus (such as systemic lupus erythematosus) and Guillain-Barr syndrome, Atopic dermatitis, autoimmune hepatitis, fibrosing alveolitis, Grave's disease, IgA nephropathy, idiopathic thrombocytopenic purpura, Meniere's disease, pemphigus, primary biliary cirrhosis, sarcoidosis, scleroderma, Wegener's granulomatosis, other autoimmune disorders, pancreatitis, trauma (surgery), graft-versus-host disease, transplant rejection, heart disease including ischaemic diseases such as myocardial infarction as well as atherosclerosis, intravascular coagulation, bone resorption, osteoporosis, osteoarthritis, periodontitis, hypochlorhydria and cancer, including breast cancer, lung cancer, gastric cancer, ovarian cancer, hepatocellular cancer, colon cancer, pancreatic cancer, esophageal cancer, head & neck cancer, kidney, and cancer, in particular renal cell carcinoma, prostate cancer, liver cancer, melanoma, sarcoma, myeloma, neuroblastoma, placental choriocarcinoma, cervical cancer, and thyroid cancer, and the metastatic forms thereof.

[0304] The present disclosure also provides a pharmaceutical or diagnostic composition comprising a bispecific protein complex of the present invention in combination with one or more of a pharmaceutically acceptable excipient, diluent or carrier. Accordingly, provided is the use of a bispecific protein complex of the invention for use in treatment and in the manufacture of a medicament.

[0305] The composition will usually be supplied as part of a sterile, pharmaceutical composition that will normally include a pharmaceutically acceptable carrier. A pharmaceutical composition of the present invention may additionally comprise a pharmaceutically-acceptable adjuvant.

[0306] The present invention also provides a process for preparation of a pharmaceutical or diagnostic composition comprising adding and mixing the antibody molecule or bispecific antibody complex of the present invention together with one or more of a pharmaceutically acceptable excipient, diluent or carrier.

[0307] The term "pharmaceutically acceptable excipient" as used herein refers to a pharmaceutically acceptable formulation carrier, solution or additive to enhance the desired characteristics of the compositions of the present disclosure. Excipients are well known in the art and include buffers (e.g., citrate buffer, phosphate buffer, acetate buffer and bicarbonate buffer), amino acids, urea, alcohols, ascorbic acid, phospholipids, proteins (e.g., serum albumin), EDTA, sodium chloride, liposomes, mannitol, sorbitol, and glycerol. Solutions or suspensions can be encapsulated in liposomes or biodegradable microspheres. The formulation will generally be provided in a substantially sterile form employing sterile manufacture processes.

[0308] This may include production and sterilization by filtration of the buffered solvent solution used for the formulation, aseptic suspension of the antibody in the sterile buffered solvent solution, and dispensing of the formulation into sterile receptacles by methods familiar to those of ordinary skill in the art.

[0309] The pharmaceutically acceptable carrier should not itself induce the production of antibodies harmful to the individual receiving the composition and should not be toxic. Suitable carriers may be large, slowly metabolised macromolecules such as proteins, polypeptides, liposomes, polysaccharides, polylactic acids, polyglycolic acids, polymeric amino acids, amino acid copolymers and inactive virus particles.

[0310] Pharmaceutically acceptable salts can be used, for example mineral acid salts, such as hydrochlorides, hydrobromides, phosphates and sulphates, or salts of organic acids, such as acetates, propionates, malonates and benzoates.

[0311] Pharmaceutically acceptable carriers in therapeutic compositions may additionally contain liquids such as water, saline, glycerol and ethanol. Such carriers enable the pharmaceutical compositions to be formulated as tablets, pills, dragees, capsules, liquids, gels, syrups, slurries and suspensions, for ingestion by the patient.

[0312] The bispecific protein complexes of the invention can be delivered dispersed in a solvent, e.g., in the form of a solution or a suspension. It can be suspended in an appropriate physiological solution, e.g., physiological saline, a pharmacologically acceptable solvent or a buffered solution. Buffered solutions known in the art may contain 0.05 mg to 0.15 mg disodium edetate, 8.0 mg to 9.0 mg NaCl, 0.15 mg to 0.25 mg polysorbate, 0.25 mg to 0.30 mg anhydrous citric acid, and 0.45 mg to 0.55 mg sodium citrate per 1 ml of water so as to achieve a pH of about 4.0 to 5.0. As mentioned *supra* a suspension can be made, for example, from lyophilised antibody.

[0313] A thorough discussion of pharmaceutically acceptable carriers is available in Remington's Pharmaceutical Sciences (Mack Publishing Company, N.J. 1991).

[0314] The bispecific antibody complex (or bispecific/multispecific antibody molecule of the present disclosure) may be the sole active ingredient in the pharmaceutical or diagnostic composition or may be accompanied by other active ingredients including other antibody ingredients, for example anti-TNF, anti-IL-1 β , anti-T cell, anti-IFN γ or anti-LPS antibodies, or non-antibody ingredients such as xanthines. Other suitable active ingredients include antibodies capable of inducing tolerance, for example, anti-CD3 or anti-CD4 antibodies.

[0315] In a further embodiment, the antibody, fragment or composition according to the disclosure is employed in combination with a further pharmaceutically active agent, for example a corticosteroid (such as fluticasone propionate) and/or a beta-2-agonist (such as salbutamol, salmeterol or formoterol) or inhibitors of cell growth and proliferation (such as rapamycin, cyclophosphamide, methotrexate) or alternatively a CD28 and/or CD40 inhibitor. In one embodiment the inhibitor is a small molecule. In another embodiment the inhibitor is an antibody specific to the target.

[0316] The pharmaceutical compositions suitably comprise a therapeutically effective amount of the bispecific antibody complex of the invention (or a bispecific/multispecific antibody molecule of the present disclosure).

[0317] The term "therapeutically effective amount" as used herein refers to an amount of a therapeutic agent needed to treat, ameliorate or prevent a targeted disease or condition, or to exhibit a detectable therapeutic or preventative effect. For any antibody, the therapeutically effective amount can be estimated initially either in cell culture assays or in animal models, usually in rodents, rabbits, dogs, pigs or primates. The animal model may also be used to determine the appropriate concentration range and route of administration. Such information can then be used to determine useful doses and routes for administration in humans.

[0318] The precise therapeutically effective amount for a human subject will depend upon the severity of the disease state, the general health of the subject, the age, weight and gender of the subject, diet, time and frequency of administration, drug combination(s), reaction sensitivities and tolerance/response to therapy. This amount can be determined by routine experimentation and is within the judgement of the clinician. Generally, a therapeutically effective amount will be from 0.01 mg/kg to 50 mg/kg, for example 0.1 mg/kg to 20 mg/kg.

[0319] Alternatively, the dose may be 1 to 500 mg per day such as 10 to 100, 200, 300 or 400 mg per day. Pharmaceutical compositions may be conveniently presented in unit dose forms containing a predetermined amount of an active agent of the invention.

[0320] Compositions may be administered individually to a patient or may be administered in combination (e.g. simultaneously, sequentially or separately) with other agents, drugs or hormones.

[0321] The dose at which the antibody molecule of the present invention is administered depends on the nature of the condition to be treated, the extent of the inflammation present and on whether the antibody molecule is being used prophylactically or to treat an existing condition.

[0322] The frequency of dose will depend on the half-life of the antibody molecule and the duration of its effect. If the antibody molecule has a short half-life (e.g. 2 to 10 hours) it may be necessary to give one or more doses per day. Alternatively, if the antibody molecule has a long half-life (e.g. 2 to 15 days) it may only be necessary to give a dosage once per day, once per week or even once every 1 or 2 months.

[0323] In the present disclosure, the pH of the final formulation is not similar to the value of the isoelectric point of the antibody or fragment, for if the pH of the formulation is 7 then a pI of from 8-9 or above may be appropriate. Whilst not wishing to be bound by theory it is thought that this may ultimately provide a final formulation with improved stability, for example the antibody or fragment remains in solution.

[0324] The pharmaceutical compositions of this invention may be administered by any number of routes including, but not limited to, oral, intravenous, intramuscular, intra-arterial, intramedullary, intrathecal, intraventricular, transdermal, transcutaneous (for example, see WO98/20734), subcutaneous, intraperitoneal, intranasal, enteral, topical, sublingual, intravaginal or rectal routes. Hyposprays may also be used to administer the pharmaceutical compositions of the invention.

[0325] Direct delivery of the compositions will generally be accomplished by injection, subcutaneously, intraperitoneally, intravenously or intramuscularly, or delivered to the interstitial space of a tissue. The compositions can also be administered into a specific tissue of interest. Dosage treatment may be a single dose schedule or a multiple dose schedule.

[0326] Where the product is for injection or infusion, it may take the form of a suspension, solution or emulsion in an oily or aqueous vehicle and it may contain formulatory agents, such as suspending, preservative, stabilising and/or dispersing agents. Alternatively, the bispecific protein complex (or bispecific/multispecific antibody molecule of the present disclosure) may be in dry form, for reconstitution before use with an appropriate sterile liquid. If the composition

is to be administered by a route using the gastrointestinal tract, the composition will need to contain agents which protect the antibody from degradation but which release the bispecific protein complex once it has been absorbed from the gastrointestinal tract.

[0327] A nebulisable formulation according to the present disclosure may be provided, for example, as single dose units (e.g., sealed plastic containers or vials) packed in foil envelopes. Each vial contains a unit dose in a volume, e.g., 2 ml, of solvent/solution buffer.

[0328] The term "variant" as used herein refers to peptide or protein that contains at least one amino acid sequence or nucleotide sequence alteration as compared to the amino acid or nucleotide sequence of the corresponding wild-type peptide or protein. A variant may comprise at least 80%, or 85%, or 90%, or 95%, or 98% or 99% sequence identity to the corresponding wild-type peptide or protein. However, it is possible for a variant to comprise less than 80% sequence identity, provided that the variant exhibits substantially similar function to its corresponding wild-type peptide or protein.

[0329] Antigens include cell surface receptors such as T cell or B cell signalling receptors, co-stimulatory molecules, checkpoint inhibitors, natural killer cell receptors, Immunoglobulin receptors, TNFR family receptors, B7 family receptors, adhesion molecules, integrins, cytokine/chemokine receptors, GPCRs, growth factor receptors, kinase receptors, tissue-specific antigens, cancer antigens, pathogen recognition receptors, complement receptors, hormone receptors or soluble molecules such as cytokines, chemokines, leukotrienes, growth factors, hormones or enzymes or ion channels, epitopes, fragments and post translationally modified forms thereof.

[0330] In one embodiment, the bispecific protein complex comprises one or two cell surface receptor specificities.

[0331] In one embodiment, the bispecific protein complex comprises one or two cytokine or chemokine specificities.

[0332] Antibodies or fragments to a pair of targets identified by the method according to the present disclosure may be incorporated into any format suitable for use as a laboratory reagent, an assay reagent or a therapeutic.

[0333] Thus in one aspect the disclosure extends to use of antibodies fragments or combinations thereof as pairs in any format, examples of which are given above.

[0334] The disclosure also extends to compositions, such as pharmaceutical compositions comprising said novel formats with the particular antigen specificity.

[0335] In a further aspect the disclosure includes use of the formats and the compositions in treatment.

[0336] In one embodiment, the bispecific protein complex of the present disclosure may be used to functionally alter the activity of the antigen or antigens of interest. For example, the bispecific protein complex may neutralize, antagonize or agonise the activity of said antigen or antigens, directly or indirectly.

[0337] The present disclosure also extends to a kit that is not part of the invention, for example comprising:

a) one or more fusion proteins (A-X) comprising a first Fab or Fab' (A) attached to a first binding partner of a binding pair (X); and

b) one or more fusion proteins (B-Y) comprising a second Fab or Fab' (B) attached to a second binding partner of the binding pair (Y), wherein the latter is specific for the first binding partner;

wherein the specific interaction (such as a binding interaction) of the two binding partners forms a heterodimer-tether which physically brings the two fusion proteins from a) and b) together to form a bispecific protein complex; wherein the fusion protein(s) is/are in a complexed or a non-complexed form;

1. wherein the bispecific protein complex formed the formula A-X:Y-B wherein:

A-X is a first fusion protein;

Y-B is a second fusion protein;

X:Y is a heterodimeric-tether;

: is a binding interaction between X and Y;

A is a first protein component of the bispecific protein complex selected from a Fab or Fab' fragment;

B is a second protein component of the bispecific protein complex selected from a Fab or Fab' fragment;

wherein X is independently selected from a scFv, or a peptide, with the proviso that

when X is a peptide Y is a scFv and when X is scFv then Y is a peptide;

wherein X or Y is a scFv specific to the peptide GCN4 represented by SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1; wherein the scFv is designated 52SR4 and is represented by SEQ ID NO:3 or amino acids 1 to 243 of SEQ ID NO:3 and wherein X or Y is a peptide GCN4 (SEQ ID NO:1 or amino acids 1 to 38 of SEQ ID NO:1).

[0338] Advantageously, the kit may comprise bispecific protein complexes of the present disclosure, or may comprise fusion proteins which are in a complexed or non-complexed form. In the former case, the bispecific protein complexes are ready for use "out of the box" which provides convenience and ease of use, whereas in the latter case, the bispecific

protein complexes can be assembled according to the user's requirements by combining different fusion proteins.

[0339] In another embodiment, the kit further comprises instructions for use.

[0340] In yet another embodiment, the kit further comprises one or more reagents for performing one or more functional assays.

[0341] In one embodiment, fusion proteins, bispecific proteins complexes, multiplexes, grids, libraries, compositions etc as described herein are for use as a laboratory reagent.

[0342] In a further aspect, there is provided a nucleotide sequence, for example a DNA sequence encoding a fusion protein and/or a bispecific protein complex as defined above.

[0343] In one embodiment, there is provided a nucleotide sequence, for example a DNA sequence encoding a bispecific protein complex according to the present disclosure.

[0344] In one embodiment there is provided a nucleotide sequence, for example a DNA sequence encoding a bispecific or multispecific antibody molecule according to the present disclosure.

[0345] The disclosure herein also extends to a vector comprising a nucleotide sequence as defined above.

[0346] The term "vector" as used herein refers to a nucleic acid molecule capable of transporting another nucleic acid to which it has been linked. An example of a vector is a "plasmid," which is a circular double stranded DNA loop into which additional DNA segments may be ligated. Another type of vector is a viral vector, wherein additional DNA segments may be ligated into the viral genome. Certain vectors are capable of autonomous replication in a host cell into which they are introduced (e.g., bacterial vectors having a bacterial origin of replication and episomal mammalian vectors). Other vectors (e.g., non-episomal mammalian vectors) can be integrated into the genome of a host cell, where they are subsequently replicated along with the host genome. In the present specification, the terms "plasmid" and "vector" may be used interchangeably as a plasmid is the most commonly used form of vector.

[0347] General methods by which the vectors may be constructed, transfection methods and culture methods are well known to those skilled in the art. In this respect, reference is made to "Current Protocols in Molecular Biology", 1999, F. M. Ausubel (ed), Wiley Interscience, New York and the Maniatis Manual produced by Cold Spring Harbor Publishing.

[0348] The term "selectable marker" as used herein refers to a protein whose expression allows one to identify cells that have been transformed or transfected with a vector containing the marker gene. A wide range of selection markers are known in the art. For example, typically the selectable marker gene confers resistance to drugs, such as G418, hygromycin or methotrexate, on a host cell into which the vector has been introduced. The selectable marker can also be a visually identifiable marker such as a fluorescent marker for example.

[0349] Examples of fluorescent markers include rhodamine, FITC, TRITC, Alexa Fluors and various conjugates thereof.

[0350] Also provided is a host cell comprising one or more cloning or expression vectors comprising one or more DNA sequences encoding an antibody of the present disclosure. Any suitable host cell/vector system may be used for expression of the DNA sequences encoding the antibody molecule of the present disclosure. Bacterial, for example *E. coli*, and other microbial systems may be used or eukaryotic, for example mammalian, host cell expression systems may also be used. Suitable mammalian host cells include CHO, myeloma or hybridoma cells.

[0351] The present disclosure also provides a process for the production of a fusion protein according to the present disclosure comprising culturing a host cell containing a vector of the present disclosure under conditions suitable for leading to expression of protein from DNA encoding the molecule of the present disclosure, and isolating the molecule.

[0352] The bispecific protein complexes of the present disclosure may be used in diagnosis/detection kits, wherein bispecific protein complexes with particular combinations of antigen specificities are used. For example, the kits may comprise bispecific antibody complexes that are specific for two antigens, both of which are present on the same cell type, and wherein a positive diagnosis can only be made if both antigens are successfully detected. By using the bispecific antibody complexes of the present disclosure rather than two separate antibodies or antibody fragments in a non-complexed form, the specificity of the detection can be greatly enhanced.

[0353] In one embodiment, the bispecific antibody complexes are fixed on a solid surface. The solid surface may for example be a chip, or an ELISA plate.

[0354] Further provided is the use of a bispecific protein complex of the present disclosure for detecting in a sample the presence of a first and a second peptide, whereby the bispecific complexes are used as detection agents.

[0355] The bispecific antibody complexes of the present disclosure may for example be conjugated to a fluorescent marker which facilitates the detection of bound antibody-antigen complexes. Such bispecific antibody complexes can be used for immunofluorescence microscopy. Alternatively, the bispecific antibody complexes may also be used for western blotting or ELISA.

[0356] In one embodiment, there is provided a process for purifying an antibody (in particular an antibody or fragment according to the invention).

[0357] In one embodiment, there is provided a process for purifying a fusion protein or bispecific protein complex according to the present disclosure comprising the steps: performing anion exchange chromatography in non-binding mode such that the impurities are retained on the column and the antibody is maintained in the unbound fraction. The step may, for example be performed at a pH about 6-8.

[0358] The process may further comprise an initial capture step employing cation exchange chromatography, performed for example at a pH of about 4 to 5.

[0359] The process may further comprise of additional chromatography step(s) to ensure product and process related impurities are appropriately resolved from the product stream.

[0360] The purification process may also comprise of one or more ultra-filtration steps, such as a concentration and diafiltration step.

[0361] "Purified form" as used *supra* is intended to refer to at least 90% purity, such as 91, 92, 93, 94, 95, 96, 97, 98, 99% w/w or more pure.

[0362] In the context of this specification "comprising" is to be interpreted as "including".

References

[0363]

1. Ribosome display efficiently selects and evolves high-affinity antibodies *in vitro* from immune libraries. Hanes J, Jermutus L, Weber-Bornhauser S, Bosshard HR, Plückthun A. (1998) Proc. Natl. Acad. Sci. U.S.A. 95, 14130-14135
2. Directed *in Vitro* Evolution and Crystallographic Analysis of a Peptide-binding Single Chain Antibody Fragment (scFv) with Low Picomolar Affinity. Zhand C, Spinelli S, Luginbuhl B, Amstutz P, Cambillau C, Pluckthun A. (2004) J. Biol. Chem. 279, 18870-18877
3. Antigen recognition by conformational selection. Berger C, Weber-Bornhauser S, Eggenberger Y, Hanes J, Pluckthun A, Bosshard H. R. (1999) F.E.B.S. Letters 450, 149-153

EXAMPLES

General methods employed in some of the Examples

[0364] General method 1: Human PBMC derived from platelet apheresis cones were banked as frozen aliquots. Prior to an assay being performed, cells were thawed, washed in DMEM (Life Technologies) and allowed to acclimatise to a 37 °C and 5% CO₂ environment.

[0365] General method 2: Fab'A-X and Fab'B-Y were incubated together for 90 minutes (in a 37°C/5%CO₂ environment) before mixing with 2.5x10⁵ PBMC in V-bottomed 96 well plates. PBMC plus bispecific (Fab'A-X and Fab'B-Y) or bivalent (e.g. Fab'A-X Fab'A'-Y) combinations were then incubated together for a further 90 minutes. After this time B cells were activated by the addition of 200nM of goat F(ab')₂ anti-human IgM (Southern Biotechnology) for 8 minutes at 37 °C. The signalling reaction was then halted by adding an equal volume of Cytotfix buffer (BD Biosciences). Plates were then left at room temperature for 15 minutes before centrifugation at 500g for 5 minutes. Excess supernatant was discarded from the cell pellet which was resuspended in flow buffer and washed once more. Cells were then resuspended in ice cold Perm Buffer III (BD Biosciences) for 30 minutes before being washed twice in flow buffer.

[0366] General method 3: Cells were activated as described in general method 2 and stained with a fluorescently labelled anti-CD20 antibody (BD Biosciences), anti-phospho Akt antibody that recognises a modified serine residue at position 473, an anti-phospho PLCg2 antibody that recognises a modified tyrosine residue at position 759 and an anti-IKB antibody that recognised total IκB. Plates were then resuspended and incubated for 1 hour at room temperature in the dark. After this time plates were washed a further two times and resuspended in 25 µl of flow buffer. Cellular expression of CD20, Akt and PLCg2 was measured using an Intellicyt HTFC™ flow cytometer

Example 1 Construction of a bispecific antibody complex of the present disclosure FabB-GCN4(7P14P):52SR4-FabA

[0367] Figures 2 and 4 show a representative bispecific antibody complex of the present disclosure. The bispecific antibody complex is composed of a first and second fusion protein.

[0368] The first fusion protein (A-X) includes a Fab fragment (Fab A (also referred to as Fab#1) with specificity to antigen 6, which is attached to X a scFv (clone 52SR4 SEQ ID NO: 3) via a peptide linker ASGGGG SEQ ID NO: 71 which is linked to the c-terminal of the CH₁ domain of the Fab fragment and the V_L domain of the scFv. The scFv itself also contains a peptide linker located in between its V_L and V_H domains.

[0369] The second fusion protein (B-Y) includes a Fab fragment (Fab B [also referred to as Fab#2] with specificity to antigen 5). However, in comparison to the first protein, the Fab fragment is attached to Y a peptide GCN4 (clone 7P14P SEQ ID NO: 1) via a peptide linker ASGGGG SEQ ID NO: 71 which is linked to the CH₁ domain of the Fab fragment.

[0370] The scFv, X, is specific for and complementary to the binding partner Y, GCN4. As a result, when the two fusion proteins are brought into contact with each other, a non-covalent binding interaction between the scFv and GCN4 peptide

occurs, thereby physically retaining the two fusion proteins in the form of a bispecific antibody complex.

[0371] The single chain antibody (scFv) 52SR4 was derived by constructing and panning a ribosome display VL-linker-VH scFv library from the spleens of mice immunized with GCN4(7P14P) (Reference 1). A further 2004 publication describes the affinity maturation of 52SR4 to a reported 5pM again using ribosome display of randomised libraries (Reference 2).

[0372] The GCN4 peptide was derived from the yeast transcription factor GCN4 by inclusion of Proline residues at positions 7 and 14, hence GCN4(7P14P) remains in a monomeric state on scFv binding as described in a 1999 publication by Berger *et al* (Reference 3).

[0373] The nucleotide sequences encoding the GCN4 peptide and the 52SR4 scFv were cloned into two separate vectors downstream of in-house heavy chain Fab expression vectors which contain CH₁ and which are already designed to receive antibody VH-regions.

[0374] VH-regions from an anti-antigen 6 antibody and an anti-antigen 5 antibody were then cloned separately into these two heavy chain vectors.

[0375] The nucleotide sequences encoding the GCN4 peptide and the 52SR4 scFv were separately cloned into a first and second vector respectively downstream of in-house light chain Fab expression vectors which contain CK and which are designed to receive antibody VL-regions.

[0376] VL-regions from an anti-antigen 6 antibody and an anti-antigen 5 antibody were cloned separately in frame with CK in an in-house light chain expression vector for co-expression with the appropriate heavy chain vector to express the Fab-scFv and Fab-peptide proteins. The vectors were then sequenced to confirm that the cloning has been successful and that the cells subsequently separately expressed Fab-scFv and Fab-peptide proteins with the V-regions from the anti-antigen 6 antibody and the anti-antigen 5 antibody respectively. Antigen 5 and 6 in Example 1 are not the antigens labelled antigen 5 and antigen 6 in later Examples with the large grid formats.

Example 2 - Flow cytometry demonstration of scFv:peptide interaction forming a non-covalent bispecific antibody that can co-engage two target antigens simultaneously

[0377] Figure 5 shows the results of a flow cytometry experiment which demonstrates the antigen specificities of two different bispecific antibody complexes formed using the scFv:peptide binding interaction.

[0378] The first bispecific antibody complex was constructed using the following two fusion proteins:

1. Anti-antigen 5 Fab-scFv (52SR4); and
2. Anti-antigen 6 Fab-peptide (GCN4)

[0379] The second bispecific antibody complex was constructed using the following two fusion proteins:

1. Anti-antigen 5 Fab-peptide (GCN4); and
2. Anti-antigen 6 Fab-scFv (52SR4)

[0380] Therefore, the two bispecific antibody complexes had the same Fab fragments and same binding partners (i.e. 52SR4 and GCN4). The difference between the two bispecific antibody complexes was in which Fab fragment is attached to which binding partner.

[0381] The control mixture which did not form a complex was made from the following fusion proteins:

1. Anti-antigen 5 Fab:GCN4; and
2. Anti-antigen 6 Fab:GCN4

[0382] To demonstrate the ability of the bispecific antibody complexes to bind to antigen 5, the complexes were incubated with Jurkat cells which express antigen 5. To demonstrate the ability of the bispecific antibody complexes to bind to antigen 6, the complexes once bound to antigen 5 on Jurkat cells were subsequently contacted with biotinylated antigen 6. The biotinylated antigen 6 was then detected using fluorescently labelled streptavidin.

[0383] The Jurkat cells were then run through a Facscalibur flow cytometer machine, wherein the fluorescently labelled Jurkat cells which can only be labelled when bound to a bispecific antibody complex, which is in turn bound to antigen 6, thereby indicating that the bispecific antibody complex is capable of binding to both antigen 5 and antigen 6 can be separated from Jurkat cells incubated with two fusion proteins capable of binding to antigen 5 and antigen 6, both fused to peptide which cannot form a complex.

[0384] The FACS plot in Figure 5 shows significant shifts for both the bispecific antibody complexes (thin and thick line over and above background filled), thus demonstrating that the bispecific antibody complexes can successfully bind to both target antigens and that the ability to bind to both target antigens is retained regardless of whether a given Fab

fragment is connected to a scFv or peptide.

[0385] The subsequent capture of either peptide or scFv respectively C-terminally fused to the anti-antigen 6 Fab allows further capture of biotinylated antigen 6 which is detected in a final layer with fluorescently labelled streptavidin. Accordingly, the results depicted in the FACS plot shows that the bispecific antibody complexes of the present disclosure

are able to successfully bind two different target antigens simultaneously.

[0386] Antigen 5 and 6 in Example 2 are not the antigens labelled antigen 5 and antigen 6 in later Examples with the large grid formats below.

Example 3-Biacore demonstration of scFv:peptide interaction

[0387] Figure 6 shows a surface plasmon resonance trace which demonstrates the affinity of the scFv:peptide (i.e. 52SR4:GCN4) interaction. Surface plasmon resonance was performed using a Biacore 3000 (GE Healthcare). All experiments were performed at 25°C. Streptavidin (produced in-house) was immobilised on a CM5 Sensor Chip (GE Healthcare) via amine coupling chemistry to a final level of approximately 1750 response units. HBS-N buffer (10mM HEPES pH 7.4, 0.15M NaCl; GE Healthcare) was used as the running buffer for immobilisation and peptide capture. A 5µl injection of Biotin-GCN4 peptide in HBS-N (10nM, M.W. 4360) was used to achieve approximately 6RU of capture on the immobilised streptavidin. The running buffer was switched to HBS-EP+ buffer (10mM HEPES pH 7.4, 0.15M NaCl, 3mM EDTA, 0.05% (v/v) surfactant P20; GE Healthcare) for measuring anti-GCN4 (52SR4) scFv binding kinetics. Three-fold serial dilutions of Fab-scFv (generated in-house) from 30nM, or HBS-EP+ buffer control, were injected over the immobilised GCN4 peptide (3min association, 15min dissociation) at a flow rate of 30µl/min. The surface was regenerated after each injection at a flow-rate of 10µl/min by two serial 60sec injection of 2M Guanidine-HCl. Double referenced background subtracted binding curves were analysed using the 3000 BIAeval software (version 4.1) following standard procedures. Kinetic parameters were determined from fitting the 1:1 binding model algorithm. The data demonstrate that the scFv has an affinity for the peptide of 516pM.

Example 4-Production of Fab-A (Fab-scFv [A-X]) and Fab-B (Fab-peptide [B-Y]) for functional assays

[0388] Cloning strategy: Antibody variable region DNA was generated by PCR or gene synthesis flanking restriction enzyme sites DNA sequence. These sites were *HindIII* and *XhoI* for variable heavy chains and *HindIII* and *BsiWI* for variable light chains. This makes the heavy variable region amenable to ligating into the two heavy chain vectors (pNAFH with FabB-Y and pNAFH with FabA-X) as they have complementary restriction sites. This ligates the variable region upstream (or 5') to the murine constant regions and peptide Y (GCN4) or scFv X (52SR4) creating a whole reading frame. The light chains were cloned into standard in house murine constant kappa vectors (pMmCK or pMmCK S171C) which again use the same complimentary restriction sites. The pMmCK S171C vector is used if the variable region is isolated from a rabbit. The cloning events were confirmed by sequencing using primers which flank the whole open reading frame.

[0389] Cultivating CHOSXE: Suspension CHOSXE cells were pre-adapted to CDCHO media (Invitrogen) supplemented with 2 mM (100x) glutamx. Cells were maintained in logarithmic growth phase agitated at 140 rpm on a shaker incubator (Kuner AG, Birsfelden, Switzerland) and cultured at 37°C supplemented with 8% CO₂.

[0390] Electroporation Transfection: Prior to transfection, the cell numbers and viability were determined using CEDEX cell counter (Innovatis AG, Bielefeld, Germany) and required amount of cells (2x10⁸ cells/ml) were transferred into centrifuge conical tubes and were spun at 1400 rpm for 10 minutes. The pelleted cells were re-suspended in sterile Earls Balanced Salts Solution and spun at 1400 rpm for further 10 minutes. Supernatant was discarded and pellets were re-suspended to desired cell density.

[0391] Vector DNA at a final concentration of 400 µg for 2x10⁸ cells/ml mix and 800 µl was pipetted into cuvettes (Biorad) and electroporated using in-house electroporation system. **Fab-A (Fab-scFv [A-X]) and Fab-B (Fab-peptide [B-Y]) were transfected separately.**

[0392] Transfected cells were transferred directly into 1x3L Erlenmeyer Flasks contained ProCHO 5 media enriched with 2 mM glutamx and antibiotic antimetabolic (100X) solution (1 in 500) and cells were cultured in Kuhner shaker incubator set at 37°C, 5% CO₂ and 140 rpm shaking. Feed supplement 2g/L ASF (AJINOMOTO) was added at 24hr posttransfection and temperature dropped to 32°C for further 13 days culture. At day four 3 mM sodium butyrate (n-butyric acid sodium Salt, Sigma B-5887) was added to the culture.

[0393] On day 14, cultures were transferred to tubes and supernatant separated from the cells after centrifugation for 30 minutes at 4000rpm. Retained supernatants were further filtered through 0.22µm SARTOBRA® P Millipore followed by 0.22 µm Gamma gold filters. Final expression levels were determined by Protein G-HPLC.

[0394] Large Scale (1.0L) Purification: The Fab-A and Fab-B were purified by affinity capture using the AKTA Xpress systems and HisTrap Excel pre-packed nickel columns (GE Healthcare). The culture supernatants were 0.22 µm sterile filtered and pH adjusted to neutral, if necessary, with weak acid or base before loading onto the columns. A secondary

wash step, containing 15-25 mM imidazole, was used to displace any weakly bound host cell proteins/non-specific His binders from the nickel resin. Elution was performed with 10 mM sodium phosphate, pH7.4 + 1M NaCl + 250 mM imidazole and 2 ml fractions collected. One column volume into the elution the system was paused for 10 minutes to tighten the elution peak, and consequently decrease the total elution volume. The cleanest fractions were pooled and buffer exchanged into PBS (Sigma), pH7.4 and 0.22 μ m filtered. Final pools were assayed by A280 Scan, SE-HPLC (G3000 method), SDS-PAGE (reduced & non-reduced) and for endotoxin using the PTS Endosafe system.

Example 5-Use of Fab-A (Fab-scFv [A-X]) and Fab-B (Fab-peptide [B-Y]) in heterodimerically-tethered bispecific protein complex format to select functional, bivalent and bispecific antigen target combinations based on inhibition of Akt signalling (as a measure of B cell activation)

[0395] Human PBMC were prepared according to general method 1. During this period grids of bispecific or bivalent antibodies were created by diluting equimolar (200 nM) quantities of Fab'-A (Fab-scFv) and Fab'-B (Fab-peptide) with varying antigen specificity for the cell surface proteins antigen 3, antigen 1, antigen 4 and antigen 2 in DMEM containing 10% calf serum and 2mM glutamine. This grid is shown in Table 4.

Table 4: Possible grid of bispecific and bivalent combinations of antibodies with specificity for antigen 3, antigen 1, antigen 4 and 2.

(A-X)	(B-Y) Fab B			
Fab A	Antigen 3-Y	Antigen 1-Y	Antigen 4-Y	Antigen 2-Y
Antigen 3-X	3-X:Y-3	3-X:Y-1	3-X:Y-4	3-X:Y-2
Antigen 1-X	1-X:Y-3	1-X:Y-1	1-X:Y-4	1-X:Y-2
Antigen 4-X	4-X:Y-3	4-X:Y-1	4-X:Y-4	4-X:Y-2
Antigen 2-X	2-X:Y-3	2-X:Y-1	2-X:Y-4	2-X:Y-2

where X is a scFv (52SR4) and Y is a peptide (GCN4)
 Fab'A-X and Fab'B-Y were incubated with PMBCs according to general method 2 following purification as described in Example 4.

[0396] Cells were then stained with a fluorescently labelled anti-CD20 antibody (BD Biosciences) and a fluorescently labelled anti-phospho Akt antibody that recognises a modified serine residue at position 473 on the protein. Plates were then resuspended and incubated for 1 hour at room temperature in the dark. After this time plates were washed a further two times and resuspended in 25 μ l of flow buffer. Cellular expression of CD20 and Akt was measured using an Intellicyt HTFC™ flow cytometer.

[0397] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of Akt levels was calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only). The relative effect of the combinations of antibodies to antigen 3 (VR0982), antibodies to antigen 1 (VR4247), antibodies to antigen 4 (VR4248) and antibodies to antigen 2 (VR4246) is shown in table 5 ($\downarrow$ = inhibition, $\uparrow$ =stimulation and $\leftrightarrow$ = no overall effect). The number of arrows is indicative of the intensity of the activity.

Table 5: Table of the relative potency of inhibition of phosphorylated Akt for bispecific & bivalent combinations of antibodies with specificity for antigen 3, 1, 4 & 2

(A-X) Fab A specific to	(B-Y) Fab B			
	Antigen 3-Y	Antigen 1-Y	Antigen 4-Y	Antigen 2-Y
Antigen 3-X	$\uparrow\uparrow$	Not Tested	$\uparrow\uparrow$	$\downarrow\downarrow\downarrow$
Antigen 1-X	$\leftrightarrow$	$\leftrightarrow$	$\leftrightarrow$	$\downarrow\downarrow\downarrow\downarrow$
Antigen 4-X	Not Tested	Not Tested	Not Tested	Not Tested
Antigen 2-X	$\downarrow\downarrow\downarrow$	$\downarrow\downarrow\downarrow$	$\downarrow\downarrow$	$\leftrightarrow$

where X is a scFv (52SR4) and Y is a peptide (GCN4)

[0398] This data is also shown in the form of a histogram (Figure 7): as mean values and the error bars show 95% confidence intervals. The data shows that the combinations of Fab to antigen 3 (VR0982) with Fab to antigen 2 (VR4246), Fab to antigen 1 (VR4247) with Fab to antigen 2 (VR4246) and Fab to antigen 4 (VR4248) with Fab to antigen 2 (VR4246) can all inhibit phospho-Akt expression in B cells stimulated with anti-IgM. In contrast, the combinations of Fab to antigen 3 (VR0982) with Fab to antigen 3 (VR0982) and Fab to antigen 3 (VR0982) with Fab to antigen 4 (VR4248) exhibited elevated levels of phospho-Akt expression. All other combinations tested showed no effect.

Example 6-Use of the heterodimerically-tethered bispecific protein complex format to select functional, bivalent and bispecific antigen target combinations based on inhibition of PLCg2 signalling (as a measure of B cell activation).

[0399] Human PBMC were prepared according to general method 1. During this period grids of bispecific or bivalent antibodies were created by diluting equimolar (200 nM) quantities of Fab'-A (Fab-scFv [A-X]) and Fab'-B (Fab-peptide [B-Y]) with antigen specificity for the cell surface proteins antigen 3, antigen 1, antigen 4 and antigen 2 in DMEM containing 10% calf serum and 2 mM glutamine. This grid is shown in Table 6.

[0400] Fab'A-X and Fab'B-Y were incubated according to general method 2.

[0401] Cells were then stained with a fluorescently labelled anti-CD20 antibody (BD Biosciences) and a fluorescently labelled anti-phospho PLCg2 antibody that recognises a modified tyrosine residue at position 759 on the protein. Plates were then resuspended and incubated for 1 hour at room temperature in the dark. After this time plates were washed twice and resuspended in 25 µl of flow buffer. Cellular expression of CD20 and PLCg2 was measured using an Intellicyt HTFC™ flow cytometer.

[0402] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of PLCg2 levels was calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only). The relative effect of the "antibody" combinations to antigen 3, antigen 1, antigen 4 and antigen 2 is shown in table 6 (↓= inhibition, ↑=stimulation and ↔ = no overall effect). The number of arrows is indicative of the intensity of the activity.

Table 6: Table of the relative potency of inhibition of phosphorylated PLCg2 for bispecific & bivalent combinations of antibodies with specificity for antigen 3, 1, 4 & 2.

(A-X) Fab A specific to	(B-Y) Fab B			
	Antigen 3-Y	Antigen 1-Y	Antigen 4-Y	Antigen 2-Y
Antigen 3-X	↑	Not Tested	↑↑	↓↓↓
1-X	↔	↔	↔	↓↓
4-X	Not Tested	Not Tested	Not Tested	Not Tested
2-X	↓↓↓	↓↓	↓↓↓	↔
where X is a scFv and Y is a peptide				

[0403] This data is also represented as a histogram (Figure 8), showing mean values and the error bars are 95% confidence intervals. The data shows that the combinations of a Fab to antigen 3 (VR0982) with a Fab to antigen 2 (VR4246), a Fab to antigen 1 (VR4247) with a Fab to antigen 2 (VR4246), a Fab to antigen 4 (VR4248) with a Fab to antigen 2 (VR4246) can all inhibit phospho-PLCg2 expression in B cells stimulated with anti-IgM. In contrast, the combinations of a Fab to antigen 3 (VR0982) with a Fab to antigen 3 (VR0982), and a Fab to antigen 3 (VR0982) with a Fab to antigen 4 (VR4248) exhibited elevated levels of phospho-PLCg2 expression. The combination of a Fab to antigen 1 with a Fab to antigen 1 showed no effect.

Example 7-The use of the heterodimerically-tethered bispecific protein complex format to select functional, bivalent and bispecific antigen target combinations based on inhibition of CD86 expression (as a measure of B cell activation).

[0404] Human PBMC were prepared according to general method 1. During this period grids of bispecific or bivalent antibodies were created by diluting equimolar (200 nM) quantities of Fab'-X (Fab-scFv) and Fab'-Y (Fab-peptide) with antigen specificity for the cell surface proteins antigen 3, antigen 1, antigen 4 and antigen 2 in DMEM containing 10% calf serum and 2mM glutamine. This grid is shown in Table 7.

[0405] Fab'A-X and Fab'B-Y were incubated together for 90 minutes (in a 37°C & 5%CO₂ environment) before mixing

with 2.5×10^5 PBMC in V-bottomed 96 well plates. PBMC plus bispecific or bivalent combinations were then incubated together for a further 90 minutes. After this time B cells were activated by the addition of 200 nM of goat F(ab')₂ anti-human IgM (Southern Biotechnology) for 24 hours at 37°C. After this time plates were placed on ice and washed once in ice cold flow buffer (PBS + 1% BSA + 0.01% NaN₃). Cells were then stained with a fluorescently labelled anti-CD19 antibody (BD Biosciences) and a fluorescently labelled anti-CD86 antibody and incubated on ice for 1 hour in the dark. After this time plates were washed a further two times and resuspended in 25 µl of flow buffer. Cellular expression of CD19 and CD86 was measured using an Intellicyt HTFC™ flow cytometer.

[0406] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of CD86 levels was calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only). The relative effect of the combinations of Fab to antigen 3 (VR0982), Fab to antigen 1 (VR4247), Fab to antigen 4 (VR4248) and Fab to antigen 2 (VR4246) is shown in table 7 (↓= inhibition, ↑=stimulation and ↔ = no overall effect). The number of arrows is indicative of the intensity of the activity.

Table 7: Table of the relative potency of inhibition of B Cell CD86 expression for bispecific and bivalent combinations of antibodies with specificity for antigen 3, 1, 4 & 2

(A-X) Fab A specific to	(B-Y) Fab B			
	Antigen 3-Y	Antigen 1-Y	Antigen 4-Y	Antigen 2-Y
Antigen 3-X	↑	Not Tested	↔	↓↓↓
Antigen 1-X	↔	↔	↑↑	↔
Antigen 4-X	Not Tested	Not Tested	Not Tested	Not Tested
Antigen 2-X	↓↓↓	↓	↓↓↓↓	↓↓

where X is a scFv (52SR4) and Y is a peptide (GCN4)

[0407] This data is also shown as a histogram (Figure 9), as mean values and the error bars are 95% confidence intervals. The data shows that the combinations of a Fab to antigen 3 (VR0982) with a Fab to antigen 2 (VR4246), a Fab to antigen 1 (VR4247) with a Fab to antigen 2 (VR4246), a Fab to antigen 4 (VR4248) with a Fab to antigen 2 (VR4246) and a Fab to antigen 2 (VR4246) with a Fab to antigen 2 (VR4246) can all inhibit CD86 expression on B cells stimulated with anti-IgM. In contrast the combinations of a Fab to antigen 3 (VR0982) with a Fab to antigen 3 (VR0982), and a Fab to antigen 1 (VR4247) with a Fab to antigen 4 (VR4248) exhibited elevated levels of CD86 expression. All the other combinations tested showed no effect.

Example 8 - The inhibitory effect of a Fab to antigen 1 (VR4247) and a Fab to antigen 2 (VR4246) can only be reproduced when the antibodies are arranged in a bispecific orientation.

[0408] Human PBMC were prepared according to general method 1. During this period combinations of bispecific, bivalent or mixtures of antibodies were created by diluting equimolar (200nM) quantities of Fab'A-X (Fab-scFv) and/or Fab'B-Y (Fab-peptide) with antigen specificity for the cell surface proteins antigen 1 and antigen 2 in DMEM containing 10% calf serum and 2mM glutamine. In addition single fab controls (Fab'-X and Fab'-Y) were also added. These combinations are shown in Table 8.

Table 8: Grid of bispecific, bivalent, mixtures or single Fab's with specificity for antigen 1 and antigen 2

(A-X) Fab A specific to	(B-Y) Fab B			
	-	Antigen 1-Y	Antigen 2-Y	Antigen 2-X
Antigen 1-X	1-X	1-X:Y-1	1-X:Y-2	1-X X-2
Antigen 2-X	2-X	2-X:Y-1	2-X:Y-2	-
Antigen 1-Y	1-Y	-	1-Y Y-2	-
Antigen 2-Y	2-Y	-	-	-

where X is a scFv (52SR4) and Y is a peptide (GCN4)
 Fab'A-X and/or Fab'B-Y were incubated according to general method 2.

[0409] Cells were then stained with a fluorescently labelled anti-CD20 antibody (BD Biosciences) and a fluorescently labelled anti-phospho Akt antibody that recognises a modified serine residue at position 473 on the protein. Plates were then resuspended and incubated for 1 hour at room temperature in the dark. After this time plates were washed a further two times and resuspended in 25 μ l of flow buffer. Cellular expression of CD20 and Akt was measured using an Intellicyt HTFC™ flow cytometer.

[0410] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of Akt levels was calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only). Figure 10 shows that only the bispecific combination of a Fab to antigen 1 (VR4247) and a Fab to antigen 2 (VR4246) but not any other combination can modulate B cell phospho -Akt levels (the data represents mean values and the error bars are 95% confidence intervals).

Example 9 The inhibitory effect of an anti-antigen 3 (VR0982) and an anti-antigen 2 (VR4246) can only be reproduced when the antibodies are arranged in a bispecific orientation

[0411] Human PBMC were prepared according to general method 1. During this period combinations of bispecific, bivalent or mixtures of antibodies were created by diluting equimolar (200 nM) quantities of Fab'-X (Fab-scFv) and/or Fab'-Y (Fab-peptide) with antigen specificity for the cell surface proteins antigen 1 and antigen 2 in DMEM containing 10% calf serum and 2 mM glutamine. These combinations are shown in Table 9.

Table 9: Grid of bispecific, bivalent or mixtures with specificity for antigen 3 & 2

(A-X) Fab A specific to	(B-Y) Fab B		
	Antigen 3-Y	Antigen 2-Y	Antigen 2-X
Antigen 3-X	3-X:Y-3	3-X:Y-2	3-X X-2
Antigen 2-X	2-X:Y-3	2-X:Y-2	-
Antigen 3-Y	-	3-Y Y-2	-

where X is a scFv (52SR4) and Y is a peptide (GCN4)
 Fab'A-X and/or Fab'B-Y were incubated according to general method 2. Cells were then stained according to general method 3.

[0412] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of Akt and PLCg2 levels were calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only). Figures 11 and 12 show that only the bispecific combination to antigen 3 and antigen 2 but not the mixtures of an anti-antigen 3 (VR0982) and an anti-antigen 2 (VR4246) antibodies inhibited phosphorylated Akt and PLCg2 expression (the data represents mean values and the error bars are 95% confidence intervals).

[0413] In order to validate the inhibition seen with the bispecific combination to antigen 3 and to antigen 2, this combination along with a mixture of anti-antigen 3 (VR0982) and anti-antigen 2 (VR4246) antibodies was titrated and inhibition of total intracellular I κ B (signalling readout) and CD86 (activation marker after 24 hours) was measured in B cells.

[0414] As can be seen in Figure 13, a combination of antigen-3-X/antigen-2-Y but not the combination of antigen-3-X/ antigen-2-X (i.e. as a simple unlinked mixture) was able to inhibit NF- κ B signal activation after anti-IgM stimulation as measured by the level of total I κ B protein. The IC₅₀, as extrapolated using a 4 parameter logistic curve fit using Graphpad Prism 6, was 7.5 nM (the data represents mean values and the error bars are standard deviations). Additionally a titration of the combination of antigen-3-X/ antigen-2-Y but not the combination of antigen-3-X/ antigen-2-X was able to inhibit anti-IgM induced CD86 expression on B cells after 24 hours (see Figure 14). The IC₅₀, as extrapolated using a 4 parameter logistic curve fit using Graphpad Prism 6, was 10.3 nM (the data represents mean values and the error bars are standard deviations).

Example 10 The inhibitory effect of an anti-antigen 4 and an anti-antigen 2 can only be reproduced when the antibodies are arranged in a bispecific orientation.

[0415] Human PBMC were prepared according to general method 1. During this period combinations of bispecific, bivalent or mixtures of antibodies were created by diluting equimolar (200nM) quantities of Fab'A-X (Fab-scFv) and/or Fab'B-Y (Fab-peptide) with antigen specificity for the cell surface proteins antigen 4 and antigen 2 in DMEM containing

10% calf serum and 2mM glutamine. These combinations are shown in Table 10.

Table 10: Grid of bispecific, bivalent or mixtures with specificity for antigen 4 & 2

(A-X) Fab A	(B-Y) Fab B	
	antigen 2-X	antigen 2-Y
antigen 4 -Y	2-X:Y-4	4-Y Y-2
antigen 2-X	-	2-X:Y-2
where X is a scFv (52SR4) and Y is a peptide (GCN4)		

[0416] Fab'A-X and/or Fab'B-Y were incubated according to general method 2. Cells were then stained according to general method 3.

[0417] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of Akt and PLCg2 levels were calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only).

[0418] Figures 15 and 16 show that only the bispecific combination of an anti-antigen 4 (VR4248) and an anti-antigen 2 (VR4246) but not the mixtures of anti-antigen 4 (VR4248) and anti-antigen 2 (VR4246) antibodies inhibited phosphorylated Akt and PLCg2 expression (the data represents mean values and the error bars are 95% confidence intervals).

[0419] In order to validate the inhibition seen with the bispecific combination of an anti-antigen 4 (VR4248) and an anti-antigen 2 (VR4246), this combination was titrated in an assay system measuring anti-IgM induced CD86 expression on B cells.

[0420] As can be seen in Figure 17, a titration of the combination of antigen 4-X/antigen 2-Y was able to inhibit anti-IgM induced CD86 expression on B cells after 24 hours. The IC₅₀, as extrapolated using a 4 parameter logistic curve fit using Graphpad Prism 6, was 4.7 nM (the data represents mean values and the error bars are standard deviations).

Example 11 Bispecific Complex Characterisation

[0421] Purification of Functional Screening reagents: The functional screening formats Fab-X (Fab-scFv-His) and Fab-Y (Fab-peptide-His) were purified as follows after standard CHO expression. Clarified cell culture supernatants were 0.22 µm sterile filtered using a 1L stericup. The pH was measured and where necessary adjusted to pH7.4. The prepared supernatants were loaded at 5ml/min onto 5ml HisTrap Nickel Excel (GE Healthcare) columns equilibrated in 10mM Sodium phosphate, 0.5 M NaCl, pH7.4. The columns were washed with 15mM imidazole, 10mM Sodium phosphate, 0.5M NaCl, pH7.4 and then eluted with 250 mM imidazole, 10 mM Sodium phosphate, 0.5M NaCl, pH7.4. The elution was followed by absorbance at 280 nm and the elution peak collected. The peak elutions were analysed by size exclusion chromatography on a TSKgel G3000SWXL; 5µm, 7.8x300mm column developed with an isocratic gradient of 0.2M phosphate, pH7.0 at 1ml/min, with detection by absorbance at 280nm. Samples of sufficient purity were concentrated to >1m/ml and diafiltered into PBS pH7.4 (Sigma Aldrich Chemicals) using Amicon Ultra-15 concentrators with a 10kDa molecular weight cut off membrane and centrifugation at 4000xg in a swing out rotor. Where product quality was not sufficient the nickel column elutions were concentrated and applied to either a XK16/60 or XK16/60 Superdex200 (GE Healthcare) column equilibrated in PBS, pH7.4 (Sigma Aldrich Chemicals). The columns were developed with an isocratic gradient of PBS, pH7.4 (Sigma Aldrich Chemicals) at 1ml/min or 2.6ml/min respectively. Fractions were collected and analysed by size exclusion chromatography on a TSKgel G3000SWXL; 5µm, 7.8x300mm column developed with an isocratic gradient of 0.2 M phosphate, pH7.0 at 1ml/min, with detection by absorbance at 280 nm. Selected fractions were pooled and concentrated to >1mg/ml using an Amicon Ultra-15 concentrator with a 10 kDa molecular weight cut off membrane and centrifugation at 4000xg in a swing out rotor.

Analysis of Bispecific Formation in Solution

Experiment 1

[0422] Purified Fab-X (VR4247) and purified Fab-Y (VR4248) were mixed in a one to one molar ratio, with a total protein concentration of 500µg/ml and incubated overnight at ambient temperature. Controls consisted of the individual parts of the mixture at the same concentration as they would be in the mixture. 100µl of the sample and each control was injected onto a TSKgel G3000SWXL; 5 µm, 7.8x300mm column developed with an isocratic gradient of 0.2 M phosphate, pH7.0 at 1 ml/min. Detection was by absorbance at 280 nm (see Figure 18).

[0423] The size exclusion chromatograms in Figure 18 show that the Fab-X (VR4247) control has a main peak of 92%

of the total peak area with a retention time of 8.610 metric minutes. The Fab-Y (VR4248) control has a main peak of 94% of the total peak area with a retention time of 10.767 metric minutes. The retention times measured for the Fab-X and Fab-Y controls were converted to apparent molecular weight of 95kDa and 35kDa respectively by using a standard curve created from the retention times of BioRad gel filtration standards (151-1901) run under the same conditions. These apparent molecular weights are consistent with the expected apparent molecular weights for Fab-scFv and Fab-peptide molecules. The main peak for the Fab-X (VR4247)/Fab-Y (VR4248) mixture has a retention time of 9.289 metric minutes. This is converted as above to an apparent molecular weight of 187kDa. This apparent molecular weight is consistent with that expected for the pairing of one Fab-X (VR4247) with one Fab-Y (VR4248). The main peak is also 84% of the total peak area suggesting that most of the Fab-X (VR4247) and Fab-Y (VR4248) have formed the 1 to 1 bispecific protein complex. The small additional shoulder and peak that elute after the main peak are consistent with the Fab-X (VR4247) and Fab-Y (VR4248) starting materials.

Experiment 2

[0424] Purified Fab-X (VR4130) and Fab-Y (VR4131) were mixed in a one to one molar ratio, with a total protein concentration of 500 µg/ml. Aliquots of this mixture were then diluted with PBS pH7.4 to a concentration of 50 µg/ml and 5 µg/ml. Controls consisting of the individual parts of the mixture at the same concentration as they would be in the 500 µg/ml mixture were also set up. All mixtures and controls were incubated overnight at ambient temperature. 100 µl of all samples and controls were injected onto a TSKgel G3000SWXL; 5 µm, 7.8x300 mm column developed with an isocratic gradient of 0.2M phosphate, pH7.0 at 1 ml/min. Detection was by absorbance at 214nm (see Figure 19, Figure 20 and Table 11).

[0425] The size exclusion chromatograms in Figure 19 show that the Fab-X (VR4130) control has a main peak of 91% of total peak area with a retention time of 8.634 metric minutes. The Fab-Y (VR4131) control has a main peak of 97% of total peak area with a retention time of 9.361 metric minutes. The retention times measured for the Fab-X and Fab-Y controls were converted to apparent molecular weights of 109kDa and 55 kDa respectively by using a standard curve created from the retention times of BioRad gel filtration standards (151-1901) run under the same conditions. These apparent molecular weights are consistent with the expected apparent molecular weights for Fab-scFv and Fab-peptide molecules. The main peak for the Fab-X (VR4130)/Fab-Y (VR4131) mixture has a retention time of 8.016 metric minutes. This was converted as above to an apparent molecular weight of 198kDa. This apparent molecular weight is consistent with that expected for the pairing of one Fab-X (VR4130) with one Fab-Y (VR4131). The main peak is also 82% of the total peak area suggesting that most of the Fab-X (VR4130) and Fab-Y (VR4131) have formed the 1 to 1 complex. The two small peaks that elute after the main peak are consistent with the Fab-X (VR4130) and Fab-Y (VR4131) starting materials.

[0426] The size exclusion chromatograms in Figure 20 are for the Fab-X (VR4130)/Fab-Y (VR4131) 1 to 1 mixtures at 500 µg/ml, 50 µg/ml and 5 µg/ml concentration. All the traces are similar with corresponding peaks between samples having similar retention times and similar relative peak heights and areas. The percentage peak area is collated in Table 11, where the % of each peak remains fairly constant upon dilution of the mixture. This indicates that the Fab-X/Fab-Y 1 to 1 complex remains as a complex at all the concentrations tested. 75% of the Fab-X and Fab-Y are present as the 1 to 1 complex even when the mixture is diluted to 5 µg/ml which is equivalent to concentration of 40 nM for the complex.

Table 11: Size exclusion peak area data for Fab-X (VR4130)/Fab-Y (VR4131) 1:1 molar ratio mixtures at 500 µg/ml, 50 µg/ml and 5 µg/ml. Peaks were detected at an absorbance of 214 nm.

Concentrations		% Peak Area		
µg/ml	nM	Fab-X (VR4130)/Fab-Y (VR4131) 1 to 1 complex	Fab-X (VR4130)	Fab-Y (VR4131)
500	4000	82%	4%	5%
50	400	81%	11%	3%
5	40	75%	21%	3%

[0427] Hence, the results of these experiments indicate that a high proportion of the Fab-X and Fab-Y fusion proteins form the desired bispecific complexes, with a minimal proportion of monomers left over and no evidence of homodimer formation.

Example 12: Grid screening of large panels of heterodimerically tethered protein complexes to identify novel bispecific antibody targets.

[0428] *Introduction:* Following the successful validation of the bispecific format and screening method in the earlier examples the screening was expanded to a larger number of antigen pairs. A panel of antibody variable (V) region pairs to 23 different antigens expressed on B cells was generated. Using the Fab-Kd-Fab [i.e. A-X:Y-B wherein A and B are Fab fragments] format a grid of heterodimerically tethered protein complexes was formed representing multiple V region combinations of each of 315 different antigen pair combinations. These combinations were screened for their ability to modulate BCR (B cell receptor) signalling in a high through-put flow cytometry assay to select novel target pairs for intervention with a bispecific antibody.

[0429] *Immunisation:* DNA encoding selected antigens was obtained by gene synthesis or commercial sources & cloned into an expression vector with a strong constitutive promoter. Plasmid DNA was then transfected into Rab-9 rabbit fibroblast cells (ATCC® CRL-1414™) using an in-house electroporation system. Twenty four hours later cells were checked for antigen expression by flow cytometry & frozen in aliquots in liquid nitrogen until use. Up to 6 antigens were immunised per rabbit by either co-expression on the same cell or making mixtures of singly or multiple transfected cells. Rabbits were immunised with 3 doses of cells.

[0430] *Antibody discovery:* B cell cultures were prepared using a method similar to that described by Zubler *et al.* (1985). Briefly, spleen or PBMC-derived B cells from immunized rabbits were cultured at a density of approximately 2000-5000 cells per well in bar-coded 96-well tissue culture plates with 200 μ l/well RPMI 1640 medium (Gibco BRL) supplemented with 10% FCS (PAA laboratories Ltd), 2% HEPES (Sigma Aldrich), 1% L-Glutamine (Gibco BRL), 1% penicillin/streptomycin solution (Gibco BRL), 0.1% β -mercaptoethanol (Gibco BRL), 3% activated splenocyte culture supernatant and gamma-irradiated mutant EL4 murine thymoma cells (5×10^4 /well) for seven days at 37°C in an atmosphere of 5% CO₂.

[0431] The presence of antigen-specific antibodies in B cell culture supernatants was determined using a homogeneous fluorescence-based binding assay using HEK293 cells co-transfected with the antigens that the rabbits were immunized with. Screening involved the transfer of 10 μ l of supernatant from barcoded 96-well tissue culture plates into barcoded 384-well black-walled assay plates containing HEK293 cells transfected with target antigen (approximately 3000 cells/well) using a Matrix Platamate liquid handler. Binding was revealed with a goat anti-rabbit IgG Fcy-specific Cy-5 conjugate (Jackson). Plates were read on an Applied Biosystems 8200 cellular detection system.

[0432] Following primary screening, positive supernatants were consolidated on 96-well bar-coded master plates using an Aviso Onyx hit-picking robot and B cells in cell culture plates frozen at -80°C. Master plates were then screened in a homogeneous fluorescence-based binding assay on HEK293 cells transfected with antigens separately and Supera-vidin™ beads (Bangs Laboratories) coated with recombinant protein as a source of antigen. This was done in order to determine the antigen specificity for each well.

[0433] To allow recovery of antibody variable region genes from a selection of wells of interest, a deconvolution step was performed to enable identification of the antigen-specific B cells in a given well that contained a heterogeneous population of B cells. This was achieved using the Fluorescent foci method (Clargo *et al.*, 2014. *Mabs* 2014 Jan 1: 6(1) 143-159; EP1570267B1). Briefly, Immunoglobulin-secreting B cells from a positive well were mixed with either HEK293 cells transfected with target antigen or streptavidin beads (New England Biolabs) coated with biotinylated target antigen and a 1:1200 final dilution of a goat anti-rabbit Fc γ fragment-specific FITC conjugate (Jackson). After static incubation at 37°C for 1 hour, antigen-specific B cells could be identified due to the presence of a fluorescent halo surrounding that B cell. A number of these individual B cell clones, identified using an Olympus microscope, were then picked with an Eppendorf micromanipulator and deposited into a PCR tube. The fluorescent foci method was also used to identify antigen-specific B cells from a heterogeneous population of B cells directly from the bone marrow of immunized rabbits.

[0434] Antibody variable region genes were recovered from single cells by reverse transcription (RT)-PCR using heavy and light chain variable region-specific primers. Two rounds of PCR were performed, with the nested secondary PCR incorporating restriction sites at the 3' and 5' ends allowing cloning of the variable region into mouse Fab-X and Fab-Y (VH) or mouse kappa (VL) mammalian expression vectors. Heavy and light chain constructs for the Fab-X and Fab-Y expression vectors were co-transfected into HEK-293 cells using Fectin 293 (Life Technologies) or Expi293 cells using Expifectamine (Life Technologies) and recombinant antibody expressed in 6-well tissue culture plates in a volume of 5ml. After 5-7 days expression, supernatants were harvested. Supernatants were tested in a homogeneous fluorescence-based binding assay on HEK293 cells transfected with antigen and Supera-vidin™ beads (Bangs Laboratories) coated with recombinant protein or antigen transfected HEK cells. This was done to confirm the specificity of the cloned antibodies.

Production of small scale Fab A-X and Fab B-Y (Small Scale (50mL) Expi293 Transfection)

[0435] The Expi293 cells were routinely sub-cultured in Expi293™ Expression Medium to a final concentration of 0.5×10^6 viable cells / mL and were incubated in an orbital shaking incubator (Multitron, Infors HT) at 120 rpm 8% CO₂ and

37°C.

[0436] On the day of transfection cell viability and concentration were measured using an automated Cell Counter (Vi-CELL, Beckman Coulter). To achieve a final cell concentration of 2.5×10^6 viable cells / mL the appropriate volume of cell suspension was added to a sterile 250 mL Erlenmeyer shake flask and brought up to the volume of 42.5 mL by adding fresh, prewarmed Expi293™ Expression Medium for each 50 mL transfection.

[0437] To prepare the lipid-DNA complexes for each transfection a total of 50 µg of heavy chain and light chain plasmid DNAs were diluted in Opti-MEM® I medium (LifeTechnologies) to a total volume of 2.5 mL and 135 µL of ExpiFectamine™ 293 Reagent (LifeTechnologies) was diluted in Opti-MEM® I medium to a total volume of 2.5 mL. All dilutions were mixed gently and incubate for no longer than 5 minutes at room temperature before each DNA solution was added to the respective diluted ExpiFectamine™ 293 Reagent to obtain a total volume of 5 mL. The DNA-ExpiFectamine™ 293 Reagent mixtures were mixed gently and incubated for 20-30 minutes at room temperature to allow the DNA-ExpiFectamine™ 293 Reagent complexes to form.

[0438] After the DNA-ExpiFectamine™ 293 reagent complex incubation was completed, the 5 mL of DNA-ExpiFectamine™ 293 Reagent complex was added to each shake flask. The shake flasks were incubated in an orbital shaking incubator (Multitron, Infors HT) at 120 rpm, 8% CO₂ and 37°C.

[0439] Approximately 16-18 hours post-transfection, 250 µL of ExpiFectamine™ 293 Transfection Enhancer 1 (LifeTechnologies) and 2.5 mL of ExpiFectamine™ 293 Transfection Enhancer 2 (LifeTechnologies) were added to each shake flask.

[0440] The cell cultures were harvested 7 days post transfection. The cells were transferred into 50 mL spin tubes (Falcon) and spun down for 30min at 4000 rpm followed by sterile filtration through a 0.22µm Stericup (Merck Millipore). The clarified and sterile filtered supernatants were stored at 4°C. Final expression levels were determined by Protein G-HPLC.

[0441] Small Scale (50 ml) Purification: Both Fab-X and Fab-Y were purified separately by affinity capture using a small scale vacuum based purification system. Briefly, the 50 ml of culture supernatants were 0.22 µm sterile filtered before 500 µL of Ni Sepharose beads (GE Healthcare) were added. The supernatant beads mixture was then tumbled for about an hour before supernatant was removed by applying vacuum. Beads were then washed with Wash 1 (50 mM Sodium Phosphate 1 M NaCl pH 6.2) and Wash 2 (0.5 M NaCl). Elution was performed with 50 mM sodium acetate, pH4.0 + 1M NaCl. The eluted fractions buffer exchanged into PBS (Sigma), pH7.4 and 0.22µm filtered. Final pools were assayed by A280 scan, SE-UPLC (BEH200 method), SDS-PAGE (reduced & non-reduced) and for endotoxin using the PTS Endosafe system.

Screening assays

[0442] Donor PBMCs were rapidly thawed using a water bath set to 37°C, and carefully transferred to a 50 ml Falcon tube. They were then diluted dropwise to 5 ml in assay media to minimise the osmotic shock. The cells were then diluted to 20 ml carefully before adding the final media diluent to make the volume 50 ml. The cells were then spun at 500 g for 5 minutes before removing the supernatant and resuspending the cells in 1 ml media. The cells were then counted and diluted to 1.66×10^6 cells/ml before dispensing 30 µl per well into a V-bottom TC plate giving a final assay concentration of 5.0×10^4 cells/well. The cell plate was then stored covered in a 37°C, 5% CO₂ incubator until they were required, giving them a minimum of 1 hour to rest.

[0443] Fab-X and Fab-Y reagents were mixed in an equimolar ratio at 5x the final assay concentration in assay media and incubated for 90 min at 37°C, 5% CO₂. Samples were prepared in a 96-well U-bottom polypropylene plate and covered during the incubation. 10 µl of 5x Fab-KD-Fab mixture was added to the appropriate test wells containing cells and mixed by shaking at 1000 rpm for 30 sec prior to being incubated for 90 min at 37 °C, 5% CO₂.

[0444] The cells were then stimulated with 10µl of anti-human IgM. The final assay concentration of stimulus varied depending on the assay panel readouts, the three antibody cocktails A, B and C (detailed below) were stimulated at a final assay concentration of either 50 µg/ml (cocktail A & C) or 25 µg/ml (cocktail B). The assay plates were then gently mixed at 1000 rpm for 30 sec prior to incubation at 37°C, 5% CO₂ for 5min (antibody cocktail A & C) or 2 min (antibody cocktail B). The assay was stopped by adding 150 µl ice-cold BD CytoFix to all wells and incubated for 15min at RT. The fixed cells were then spun at 500 g for 5min to pellet the cells and allow removal of the supernatant using a BioTek ELx405 plate washer.

[0445] The pellet was re-suspended by vortexing the plate at 2400 rpm for 30 sec. The cells were then permeabilised at 4°C by adding 100 µl ice-cold BD Cell Permeabilisation Buffer III for 30 min. The cells were then washed in 100 µl FACS buffer and spun at 500 g for 5min. Supernatant was again removed by the ELx405 before using it to rapidly dispense 200 µl FACS Buffer to wash away any residual permeabilisation buffer. Cells were again spun at 500 g and the supernatant removed by inversion. During the preceding spin step the antibody cocktail was prepared in FACS Buffer and kept shielded from the light. The cells were then re-suspended by vortexing (2400 RPM, 30sec) before 20 µl of antibody cocktail was added to all wells and the plate shaken for 30 sec at 1000 rpm. The cells were then incubated for

60 min at RT in the dark.

[0446] The cells were then washed twice in 200 μ l FACS buffer with a 500 g spin and supernatant removed after each step. Finally the cells were re-suspended by vortexing for 30 sec at 2400 rpm before adding a final 20 μ l FACS buffer. The plate(s) were then read on the Intellicyt HTFC/ iQue instrument.

FACS Buffer = PBS + 1% BSA + 0.05% NaN₃ + 2mM EDTA

Antibody Cocktail A = 1:2 CD20 PerCp-Cy5.5 (BD Biosciences) + 1:5 PLC γ 2 AF88 + 1:10 Akt AF647 + 1:50 ERK1/2 PE (diluted in FACS buffer).

Antibody Cocktail B = 1:2 CD20 PerCp-Cy5.5 (BD Biosciences) + 1:5 Syk PE + 1:5 BLNK AF647 (diluted in FACS buffer)

Antibody Cocktail C = 1:5 CD20 PerCp-Cy5.5 (Biolegend) + 1:5 PLC γ 2 AF488 + 1:10 Akt AF647 + 1:5 Syk PE (diluted in FACS buffer)

Reagent	Supplier	Catalogue number
Anti-human IgM	Southern Biotech	2022-14
CytoFix	BD Biosciences	554655
Perm Buffer III	BD Biosciences	558050
Anti Akt (pS473) AF647	BD Biosciences	561670
Anti SYK (pY348) PE	BD Biosciences	558529
Anti PLC γ 2 (pY759) AF488	BD Biosciences	558507
Anti-BLNK(pY84) AF647	BD Biosciences	558443
Anti ERK1/2 (pT202/pY204) PE	BD Biosciences	561991
Anti-human CD20 PerCp-Cy5.5	BD Biosciences	558021
Anti-human CD20 AF488	BD Biosciences	558056
Anti-human CD20 PerCp-Cy5.5	Biolegend	340508
Phosphate Buffer Saline (PBS)	Fisher Scientific	10562765
RPMI 1640	Life Technologies	31870
Foetal Calf Serum (FCS)	Life Technologies	16140
Glutamax	Life Technologies	35050
Penicillin/ Streptomycin (P/S)	Life Technologies	15070
EDTA	Sigma	03690
Sodium Azide (NaN ₃)	Sigma	S2002
Bovine Serum Albumin (BSA)	Sigma	A1470

[0447] Fab-X + Fab-Y combinations were screened with either antibody cocktail A and B or C alone. All screens were conducted on cone cells from 2 different blood donors. Data was captured and evaluated using commercially available software tools. A total of 2500 Fab-X + Fab-Y combinations were screened to 315 different antigen combinations.

Results

[0448] The percentage inhibition of the induction of phosphorylation of BCR signalling cascade proteins by each Fab-Kd-Fab [i.e. A-X:Y-B where A and B are Fab fragments] combination was calculated, in this example looking for new combinations of antigens that inhibit B cell function, the criteria for a positive combination was set as at least 30% inhibition of at least two phospho-readouts by at least one combination of V regions. According to this threshold 11 new antigen pair combinations out of 315 examined met the required criteria. This represents a 3.5% hit rate demonstrating the importance of screening large numbers of combinations to find those of desired activity.

[0449] Figures 21-23 show the data for the antigen grid cross specificities. Values are percentage inhibition (negative value for activation) of phosphorylation of Syk, PLC γ 2 & AKT respectively and represent the mean of multiple V-region combinations evaluated. 315 different antigen combinations were tested and as can be seen the effect on BCR signalling by different combinations of antibody varied significantly from strong inhibition e.g. antigen 2 on Fab-X combined with antigen 3 and 4 on Fab-Y (69.66% and 70.4% inhibition of phospho Syk Figure 21) to activation e.g. antigen 6 on X and antigen 11 on Y (minus 118.10% phospho Syk Figure 21).

[0450] Each data point representing the mean % values represented in Figures 21-23 is shown for antigen 2 on Fab-X and antigen 3 on Fab-Y in Figure 24. In this case, 23 different combinations of different antibody V regions were evaluated. The same antigen combination but in alternative orientation, i.e. antigen 2 on Fab-Y and antigen 3) on Fab-X is shown in Figure 25. In this case, 9 different combinations of different antibody V-regions were evaluated. All V regions show inhibition but advantageously this method can also be used in the selection of optimal V-region combinations.

[0451] Similarly, each data point representing the mean % values represented in Figures 21-23 is shown for antigen combination 2 on Fab-X and antigen 4 on Fab-Y in figure 26. In this case, 10 different combinations of different antibody V regions were evaluated. The same antigen combination but in alternative orientation, i.e. antigen 2 on Fab-Y and antigen 4 on Fab-X is shown in Figure 27. In this case, 6 different combinations of different antibody V regions were evaluated. Again, all V regions show inhibition but optimal V region combinations can be identified and selected using the method.

Example 13 - Evaluation of transiently expressed heterodimerically tethered protein complexes to evaluate whether FabA-X:Y-FabB grid screening can identify novel bispecific antibody targets without recourse to protein purification

[0452] *Introduction:* V-regions to 2 different antigens, 2 and 3 that inhibit B cell signalling as a bispecific antibody which were identified using the Fab-Kd-Fab [FabA-X:Y-FabB] format and grid screening of heterodimerically tethered protein complexes were expressed transiently as FabA-X and FabB-Y. The activity of transiently expressed (without subsequent purification) and a purified FabA-X and FabB-Y combination was compared to evaluate whether grid screening could be conducted with the direct products of transient expression instead of purified components.

[0453] *Immunisation:* The preparation of antigen expressing cells and immunisation of rabbits was carried out in the same way as described in Example 12.

Antibody discovery

[0454] B cell cultures were prepared in the same way as described in Example 12.

[0455] The screening of antigen-specific antibodies in B cell culture supernatants and the deconvolution step for identification of antigen specific B cells was determined in the same way as described in Example 12.

[0456] Antibody variable region genes were recovered from single cells by reverse transcription (RT)-PCR using heavy and light chain variable region-specific primers. Two rounds of PCR were performed, with the nested 2^o PCR incorporating restriction sites at the 3' and 5' ends allowing cloning of the variable region into mouse Fab-X and Fab-Y (VH) or mouse kappa (VL) mammalian expression vector. A 3^o PCR was then performed enabling the combination of amplified variable regions, human CMV promoter fragment and rabbit gamma 1 heavy constant or rabbit kappa constant fragment to generate separate heavy and light transcriptionally active PCR (TAP) fragments. These DNA fragments were used directly for recombinant expression of rabbit full-length IgG antibodies in HEK-293 cells using 293Fectin (Life Technologies) or in Expi293 cells using Expifectamine (Life Technologies).

[0457] The resulting recombinant antibodies were then screened for antigen binding using a homogeneous fluorescence-based binding assay on HEK-293 cells transfected with antigen and SuperavidinTM beads (Bangs Laboratories) coated with recombinant protein. Once specificity was confirmed with the TAP transients, antibody genes were cloned into Fab-X and Fab-Y expression vectors. Heavy and light chain constructs were co-transfected into HEK-293 cells using Fectin 293 (Life Technologies) or Expi293 cells using Expifectamine (Life Technologies) and recombinant antibody expressed in 6-well tissue culture plates in a volume of 5ml. After 5-7 days expression, supernatants were harvested.

Supernatants were tested in a homogeneous fluorescence-based binding assay on HEK293 cells transfected with antigen and Superavidin™ beads (Bangs Laboratories) coated with recombinant protein or antigen transfected HEK cells. This was done to confirm the specificity of the cloned antibodies.

5 ***Production of transient supernatants containing Fab-X and Fab-Y***

[0458] The same Expi293 transfection method described in Example 12 was used to produce the transient supernatants containing Fab-X and Fab-Y.

10 ***Production of purified Fab-X and Fab-Y***

[0459] Suspension CHOSXE cells were pre-adapted to CDCHO media (Invitrogen) supplemented with 2mM (100x) glutamx. Cells were maintained in logarithmic growth phase agitated at 140rpm on a shaker incubator (Kuner AG, Birsfelden, Switzerland) and cultured at 37°C supplemented with 8% CO₂.

15 **[0460]** Prior to transfection, the cell numbers and viability were determined using CEDEX cell counter (Innovatis AG, Bielefeld, Germany) and required amount of cells (2x10⁸ cells/ml) were transferred into centrifuge conical tubes and were spun at 1400 rpm for 10 minutes. The Pelleted cells were re-suspended in sterile Earls Balanced Salts Solution and spun at 1400 rpm for further 10 minutes. Supernatant was discarded and pellets were re-suspended to desired cell density.

20 **[0461]** Vector DNA at a final concentration of 400ug for 2x10⁸ cells/ml mix and 800 µl was pipetted into Cuvettes (Biorad) and electroporated using in-house electroporation system.

[0462] Transfected cells were transferred directly into 1X3L Erlenmeyer Flasks contained ProCHO 5 media enriched with 2 mM glutamx and antibiotic antimetabolic (100x) solution (1 in 500) and Cells were cultured in Kuhner shaker incubator set at 37°C, 5% CO₂ and 140 rpm shaking. Feed supplement 2 g/L ASF (AJINOMOTO) was added at 24 hr post transfection and temperature dropped to 37°C for further 13 days culture. At day four 3 mM sodium butyrate (n-butyric acid sodium salt, Sigma B-5887) was added to the culture.

[0463] On day 14, cultures were transferred to tubes and supernatant separated from the cells after centrifugation for 30 minutes at 4000 rpm. Retained supernatants were further filtered through 0.22 µm SARTOBRA® P Millipore followed by 0.22 µm Gamma gold filters. Final expression levels were determined by Protein G-HPLC.

30 **[0464]** The Fab-X and Fab-Y were purified by affinity capture using the AKTA Xpress systems and HisTrap Excel pre-packed nickel columns (GE Healthcare). The culture supernatants were 0.22 µm sterile filtered and pH adjusted to neutral, if necessary, with weak acid or base before loading onto the columns. A secondary wash step, containing 15-25 mM Imidazole, was used to displace any weakly bound host cell proteins / non-specific His binders from the nickel resin. Elution was performed with 10mM sodium phosphate, pH7.4 + 1 M NaCl + 250 mM imidazole and 2ml fractions collected. 35 One column volume into the elution the system was paused for 10 minutes to tighten the elution peak, and consequently decrease the total elution volume. The cleanest fractions were pooled and buffer exchanged into PBS (Sigma), pH7.4 and 0.22 µm filtered. Final pools were assayed by A280 Scan, SE-HPLC (G3000 method), SDS-PAGE (reduced & non-reduced) and for endotoxin using the PTS Endosafe system.

40 ***Functional assays***

[0465] Activation Marker Assay: Antigen 2-specific Fab'-Y and antigen 3-specific Fab'-X, either purified or in transient supernatant, were incubated together for 60 minutes (in a 37°C & 5% CO₂ environment) at equimolar concentration. The combinations were titrated from a starting molarity of 185 nM, in 1:4 serial dilutions. A mock supernatant was also 45 included, titrated from neat. In V-bottomed 96 well plates, 1.5x10⁵ PBMC were added to wells, to which were added titrated Fab'-X and Fab'-Y combinations or mock supernatant. The combinations and cells were then incubated together for a further 90 minutes. After this time B cells were activated by the addition of 12.5 µg / mL of goat F(ab')₂ anti-human IgM (Southern Biotechnology) for 24 hours at 37°C plus 5% CO₂.

[0466] To the wells were added 100 µL ice-cold FACS buffer (PBS + 1% BSA + 0.1% NaN₃ + 2 mM EDTA), the plates were sealed and covered with wet-ice for approximately 15 minutes, before centrifuging at 500xg for 5 minutes at 4°C. Excess supernatant was discarded from the cell pellets and the plates shaken at 2000 rpm for 30 seconds.

[0467] Cells were then stained with a cocktail of fluorescently labelled anti-CD19, anti-CD20 and anti-CD71 antibodies (BD Biosciences). Plates were shaken briefly and incubated for 1 hour on wet-ice in the dark. After this time plates were washed twice and resuspended in 20 µL of FACS buffer. Cellular expression of CD19, CD20 and CD71 was measured 55 using an Intellicyt iQUE® Screener flow cytometer.

[0468] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of CD71 levels were calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only).

[0469] *PhosFlow Assay:* Antigen 2-specific Fab'-Y and antigen 3-specific Fab'-X, either purified or in transient supernatant, were incubated together for 60 minutes (in a 37°C & 5% CO₂ environment) at equimolar concentration. The combinations were titrated from a starting molarity of 185 nM, in 1:4 serial dilutions. A mock supernatant was also included, titrated from neat. In V-bottomed 96 well plates, 5.0x10⁴ PBMC were added to wells, to which were added titrated Fab'-X and Fab'-Y combinations or mock supernatant. The combinations and cells were then incubated together for a further 90 minutes. After this time B cells were activated by the addition of 25 µg / mL of goat F(ab')₂ anti-human IgM (Southern Biotechnology) for 15 minutes at 37°C plus 5% CO₂. The signalling reaction was then halted by adding an equal volume of Cytotfix buffer (BD Biosciences). Plates were then left at room temperature for 15 minutes before centrifugation at 500 xg for 5 minutes. Excess supernatant was discarded from the cell pellet which was resuspended in FACS buffer (PBS + 1% BSA + 0.01% NaN₃ + 2 mM EDTA) and washed once more. Cells were then resuspended in ice cold Perm Buffer III (BD Biosciences) for 30 minutes before being washed twice in flow buffer.

[0470] Cells were then stained with a fluorescently labelled anti-CD20 antibody (BD Biosciences) and an anti-phosphorylated p38 antibody that recognises the conserved dual phosphorylated site pT180/pY182. Plates were then resuspended and incubated for 1 hour at room temperature in the dark. After this time plates were washed a further two times and resuspended in 20 µL of FACS buffer. Cellular expression of CD20 and phosho-p38 was measured using an Intellicyt IQUE® flow cytometer.

[0471] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of p38 levels were calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only).

Results

[0472] *Activation Marker Assay:* As can be seen in Figure 28, the data shows that the combination of Antigen 3 with Antigen 2 whether purified or from transient supernatant, can inhibit CD71 expression on B-cells stimulated with anti-IgM.

[0473] *PhosFlow Assay:* The data in Figure 29 shows that the combination of Antigen 3 with Antigen 2 whether purified or from transient supernatant, can inhibit phosphorylated p38 in B-cells stimulated with anti-IgM.

[0474] The surprising ability to be able to construct the bispecific complexes of the present invention directly from transiently expressed cultures without recourse to purification allows even higher throughput screening of bispecific complexes to be achieved than when purified components are used.

Example 14 - Screening of transiently expressed V-regions to Antigen 3 as Fab-X with purified anti-antigen 2 Fab-Y in heterodimerically tethered protein complexes to select optimal Antigen 3 antibody V-regions

[0475] *Introduction:* New V-regions to Antigen 3 that inhibit B cell signalling as a bispecific antibody in combination with Antigen 2 specific V regions were identified using grid screening of heterodimerically tethered protein complexes. The Antigen 3 V regions were expressed transiently as Fab-X and combined with purified anti-antigen 2 Fab-Y. The inhibition of activation of B cell signalling was measured to select the most potent Antigen 3 and antigen 2 V regions.

[0476] The preparation of antigen expressing cells and immunisation of rabbits was carried out in the same way as described in Example 12.

[0477] *Antibody discovery:* B cell cultures were prepared in the same way as described in Example 12.

[0478] The screening of antigen-specific antibodies in B cell culture supernatants and the deconvolution step for identification of antigen specific B cells was determined in the same way as Example 12.

[0479] Additional variable regions were discovered by the direct foci method directly from spleen and bone marrow-derived B cells of immunized mice. Briefly, cells at a final density of between 4x10⁵ and 8x10⁵ cells/ml were mixed with streptavidin beads (New England Biolabs) coated with biotinylated antigen and a 1:1200 final dilution of goat anti-mouse Fcγ fragment-specific FITC conjugate (Jackson). After static incubation at 37°C for 1 hour, antigen-specific B cells could be identified due to the presence of a fluorescent halo surrounding that B cell. A number of these individual B cell clones, identified using an Olympus microscope, were then picked with an Eppendorf micromanipulator and deposited into a PCR tube.

[0480] Antibody variable region genes were recovered from single cells by reverse transcription (RT)-PCR using heavy and light chain variable region-specific primers. Two rounds of PCR were performed, with the nested 2nd PCR incorporating restriction sites at the 3' and 5' ends allowing cloning of the variable region into mouse Fab-X and mouse kappa (VL) mammalian expression vector. These vectors were then co-transfected in HEK-293 cells using 293Fectin (Life Technologies) or in Expi293 cells using Expifectamine (Life Technologies) and left to express for 6 days. Supernatants were tested in a homogeneous fluorescence-based binding assay on HEK293 cells transfected with antigen and Superavidin™ beads (Bangs Laboratories) coated with recombinant protein or antigen transfected HEK cells. This was done to confirm the specificity of the cloned antibodies.

[0481] In addition to the Fab-X transient supernatants, negative control Mock supernatants were prepared in the same

way using an irrelevant control DNA.

[0482] The expression levels of Fab-X were determined by Protein G-HPLC.

[0483] Production of purified Fab-Y: Purified Fab-Y was prepared using the same method described in Example 13

5 Functional assay

[0484] The same functional assay as described in Example 12 was used, except that instead of 3 different antibody cocktails, only one cocktail was used with the same assay concentrations and incubation conditions as described for antibody cocktail A in Example 12.

Antibody Cocktail = 1:3 CD20 PerCp-Cy5.5 + 1:5 PLCγ2 AF88 + 1:10 Akt AF647 + 1:5 p38 MAPK PE (diluted in FACS buffer).

15 Results

[0485] As can be seen in Figures 30-33, the data shows that the combination of different transiently expressed antigen 3 mouse V regions in Fab-X with 2 different purified antigen 2 V regions (VR447 and VR4450) in Fab-Y can inhibit B cell activation to different levels and screening therefore facilitates selection of optimal V regions. Combinations with transient Fab-X are compared to a reference combination with a purified Fab-X (VR4126).

25 Example 15 - Comparison of the activity of antigen 2 plus antigen 3 co-targeting in Fab-Kd-Fab screening format to a molecularly linked bispecific BYbe format.

[0486] Introduction: To check that target pair activity identified in the Fab-Kd-Fab heterodimerically tethered screening complex could translate to similar desired activity in an alternative therapeutic molecularly linked format, Antigen 2 specificity (VR4447) and antigen 3 specificity (VR4130) were generated in a BYbe format. This BYbe format consists of the anti-Antigen 3 V regions (VR4130) as a disulphide stabilised (ds) single chain (sc)-Fv fused to the heavy chain of the anti-Antigen 2 Fab (VR4447).

30 Methods:

[0487] As described in Example 13 except the purification of BYbes for functional screening was performed as follows: The functional screening BYbe (Fab-dsscFv [scFv off C-terminus of Fab heavy chain]) formats were purified as follows. Clarified cell culture supernatants from standard expiHEK or CHO expression were 0.22µm sterile filtered. The filtered supernatants were loaded at 2ml/min onto 50ml GammabindPlus Sepharose XK26 columns (GE Healthcare) equilibrated in PBS pH7.4 (Sigma Aldrich Chemicals). After loading the columns were washed with PBS pH7.4 and then eluted with 0.1M Glycine/HCl, pH2.7. The elution was followed by absorbance at 280nm, the elution peak collected, and then neutralised with 1/25th volume of 2M Tris/HCl pH8.5. The neutralised samples were concentrated using Amicon Ultra-15 concentrators with a 10kDa (BYbes) molecular weight cut off membrane and centrifugation at 4000xg in a swing out rotor. Concentrated samples were applied to either a XK16/60 or XK26/60 Superdex200 column (GE Healthcare) equilibrated in PBS, pH7.4. The columns were developed with an isocratic gradient of PBS, pH7.4 at either 1ml/min or 2.6ml/min respectively. Fractions were collected and analysed by size exclusion chromatography on a TSK gel G3000SWXL; 5µm, 7.8 X 300 mm column developed with an isocratic gradient of 0.2M phosphate, pH7.0 at 1 ml/min, with detection by absorbance at 280 nm. Selected monomer fractions were pooled and concentrated to >1 mg/ml using an Amicon Ultra-15 concentrator with a 10kDa molecular weight cut off membrane and centrifugation at 4000xg in a swing out rotor. Final samples were assayed; for concentration by A280 Scanning UV-visible spectrophotometer (Cary 50Bio); for % monomer by size exclusion chromatography on a TSK gel G3000SWXL; 5 µm, 7.8x300 mm column developed with an isocratic gradient of 0.2 M phosphate, pH7.0 at 1ml/min, with detection by absorbance at 280nm; by reducing and non-reducing SDS-PAGE run on 4-20% Tris-Glycine 1.5 mm gels (Novex) at 50 mA (per gel) for 53minutes; and for endotoxin by Charles River's EndoSafe® Portable Test System with Limulus Amebocyte Lysate (LAL) test cartridges.

55 Functional assays

[0488] Activation Marker Assay: Antigen 2-specific Fab'-Y and Antigen 3-specific Fab'-X, were incubated together for 60 minutes (in a 37°C and 5% CO₂ environment) at equimolar concentration. The combinations were titrated from a

starting molarity of 100 nM, in 1:4 serial dilutions. Antigen 2 and 3-specific BYbe was also titrated from a starting molarity of 100 nM, in 1:4 serial dilutions. In V-bottomed 96 well plates, 1.5×10^5 PBMC were added to wells, to which were added titrated Fab'-X and Fab'-Y combinations or titrated BYbe. The Fab'-X and Fab'-Y combinations or BYbe were incubated with cells for a further 90 minutes. After this time B cells were activated by the addition of 25 μ g / mL of goat F(ab')₂

anti-human IgM (Southern Biotechnology) for 24 hours at 37°C plus 5% CO₂.
[0489] To the wells were added 100 μ L ice-cold FACS buffer (PBS + 1% BSA + 0.1% NaN₃ + 2 mM EDTA), the plates were sealed and covered with wet-ice for approximately 15 minutes, before centrifuging at 500 xg for 5 minutes at 4°C. Excess supernatant was discarded from the cell pellets and the plates shaken at 2000 rpm for 30 seconds.

[0490] Cells were then stained with a cocktail of fluorescently labelled anti-CD19, anti-CD20 and anti-CD71, anti-CD40 and anti-CD86 antibodies (BD Biosciences). Plates were shaken briefly and incubated for 1 hour on wet-ice in the dark. After this time plates were washed twice and resuspended in 20 μ L of FACS buffer. Cellular expression of CD19, CD20 and CD71, CD40 and CD86 was measured using an Intellicyt iQUE® Screener flow cytometer. Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of CD71, CD40 and CD86 levels were calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only).

[0491] PhosFlow Assay: Antigen 2-specific Fab'-Y and Antigen 3-specific Fab'-X, were incubated together for 60 minutes (in a 37°C and 5% CO₂ environment) at equimolar concentration. The combinations were titrated from a starting molarity of 100 nM, in 1:4 serial dilutions. Antigen 2 and Antigen 3-specific BYbe was also titrated from a starting molarity of 100 nM, in 1:4 serial dilutions. In V-bottomed 96 well plates, 5.0×10^4 PBMC were added to wells, to which were added titrated Fab'-X and Fab'-Y combinations or titrated BYbe. The Fab'-X and Fab'-Y combinations or BYbe were incubated with cells for a further 90 minutes. After this time B cells were activated by the addition of 25 μ g / mL of goat F(ab')₂ anti-human IgM (Southern Biotechnology) for 15 minutes at 37°C plus 5% CO₂. The signalling reaction was then halted by adding an equal volume of Cytofix buffer (BD Biosciences). Plates were then left at room temperature for 15 minutes before centrifugation at 500 x g for 5 minutes. Excess supernatant was discarded from the cell pellet which was resuspended in FACS buffer (PBS + 1% BSA + 0.01% NaN₃ + 2 mM EDTA) and washed once more. Cells were then resuspended in ice cold Perm Buffer III (BD Biosciences) for 30 minutes before being washed twice in flow buffer.

[0492] Cells were then stained with a fluorescently labelled anti-CD20 antibody (BD Biosciences) and anti-phosphorylated PLCy2, anti-phosphorylated Akt and anti-phosphorylated p38 antibodies (BD Biosciences). Plates were then resuspended and incubated for 1 hour at room temperature in the dark. After this time plates were washed a further two times and resuspended in 20 μ L of FACS buffer. Cellular expression of CD20 and phospho-PLCy2, phospho-Akt and phospho-p38 were measured using an Intellicyt iQUE® flow cytometer.

[0493] Using the data analysis software package Forecyt™ (Intellicyt) B cells were identified as distinct from other cell populations and the geometric mean of PLCy2, Akt and p38 levels were calculated for each well. All data was then expressed as the percentage inhibition of the maximal response (anti-IgM only) minus the background (cells only).

Results

[0494] PhosFlow Assay: The data in Figure 34 show that targeting antigen 3 and antigen 2 either in the Fab-Kd-Fab or BYbe format can inhibit phosphorylated PLCy2 in B-cells stimulated with anti-IgM. The data in Figure 35 show that targeting antigen 3 and antigen 2 either in the Fab-Kd-Fab or BYbe format can inhibit phosphorylated P38 in B-cells stimulated with anti-IgM. The data in Figure 36 show that targeting antigen 3 and antigen 2 either in the Fab-Kd-Fab or BYbe format can inhibit phosphorylated Akt in B-cells stimulated with anti-IgM.

[0495] Activation Marker Assay: As can be seen in Figure 37, the data show that targeting antigen 3 and antigen 2 either in the Fab-Kd-Fab or BYbe format can inhibit CD71 expression on B-cells stimulated with anti-IgM. The data in Figure 38 show that targeting antigen 3 and antigen 2 either in the Fab-Kd-Fab or BYbe format can inhibit CD40 expression on B-cells stimulated with anti-IgM. The data in Figure 39 show that targeting antigen 3 and antigen 2 either in the Fab-Kd-Fab or BYbe format can inhibit CD86 expression on B-cells stimulated with anti-IgM.

Example 16 - Comparison of the activity of antigen 2 plus antigen 3 co-targeting in a molecularly linked bispecific Bybe format with the further addition of an anti-albumin binding domain for extension of *in vivo* half-life.

[0496] Introduction: To check that target pair activity identified in the Fab-Kd-Fab heterodimerically tethered screening complex could translate to similar desired activity in a potential therapeutic molecularly linked format with an anti-albumin targeted *in vivo* half-life extension, an anti-albumin antibody fragment was fused to the light chain of the antigen 3 Fab of the BYbe format described in Example 15. Antigen 2 specificity (VR4447) and antigen 3 specificity (VR4130 and VR4126) were generated in a Bybe format with and without addition of an anti-albumin fragment (VR0645).

[0497] Description of constructs used in this experiment.

Construct Name	Fab Specificity	Heavy Chain	Light Chain
VR4447/VR4126 BYbe	Antigen 2	Antigen 3	None
VR4447/VR4126/VR645) BYbe/ Albumin	Antigen 2	Antigen 3	Albumin
VR4447/VR4130 BYbe	Antigen 2	Antigen 3	None
VR4447/VR4130/VR645) BYbe/ Albumin	Antigen 2	Antigen 3	Albumin

Methods

Purification of BYbes for Functional Screening:

[0498] The functional screening BYbe (Fab-dsscFv [scFv off C-terminus of Fab heavy chain]) format was purified as follows. Clarified cell culture supernatants from standard expiHEK or CHO expression were 0.22 μ m sterile filtered. The filtered supernatants were loaded at 2 ml/min onto 50 ml GammabindPlus Sepharose XK26 columns (GE Healthcare) equilibrated in PBS pH7.4 (Sigma Aldrich Chemicals). After loading the columns were washed with PBS pH7.4 and then eluted with 0.1M Glycine/HCl, pH 2.7. The elution was followed by absorbance at 280nm, the elution peak collected, and then neutralised with 1/25th volume of 2 M Tris/HCl pH8.5. The neutralised samples were concentrated using Amicon Ultra-15 concentrators with either a 10 kDa or 30 kDa molecular weight cut off membrane and centrifugation at 4000 xg in a swing out rotor. Concentrated samples were applied to either a XK16/60 or XK26/60 Superdex 200 column (GE Healthcare) equilibrated in PBS, pH7.4. The columns were developed with an isocratic gradient of PBS, pH7.4 at either 1 ml/min or 2.6 ml/min respectively. Fractions were collected and analysed by size exclusion chromatography on a TSK gel G3000SWXL; 5 μ m, 7.8 X 300mm column developed with an isocratic gradient of 0.2 M phosphate, pH 7.0 at 1 ml/min, with detection by absorbance at 280 nm. Selected monomer fractions were pooled and concentrated to >1 mg/ml using an Amicon Ultra-15 concentrator with a 10 kDa or 30 kDa molecular weight cut off membrane and centrifugation at 4000 xg in a swing out rotor. Final samples were assayed; for concentration by A280 Scanning UV-visible spectrophotometer (Cary 50Bio); for % monomer by size exclusion chromatography on a TSK gel G3000SWXL; 5 μ m, 7.8x300 mm column developed with an isocratic gradient of 0.2 M phosphate, pH7.0 at 1 ml/min, with detection by absorbance at 280 nm; by reducing and non-reducing SDS-PAGE run on 4-20% Tris-Glycine 1.5 mm gels (Novex) at 50 mA (per gel) for 53 minutes; and for endotoxin by Charles River's EndoSafe® Portable Test System with Limulus Amebocyte Lysate (LAL) test cartridges.

[0499] 100 nM of each construct purified protein were pre-incubated with human PBMC derived from five separate donors for 60 min at 37 degree C/5%CO₂ in RPMI 1640 media plus 10% foetal bovine serum and 2 mM Glutamax (R10 media). After 60 min cells were stimulated with 25 ug/ml of a goat anti-IgM antibody designed to stimulate B cells only. 24 hours later plates were placed on ice to halt any further cell activation before washing once with ice cold flow cytometry buffer (PBS+1%BSA+0.01%NaN₃). All supernatant was removed and cell pellets resuspended. Cells were placed on ice and a cocktail of anti-CD19, -CD20, -CD27, -CD71 and CD86 antibodies added. Cells were incubated for 60 min before washing twice in flow cytometry buffer. Data on the binding of anti-CD27, -CD71 and CD86 to CD19/CD20 positive B cells was generated using an iQUE high throughput flow cytometer. Forecyt software was used to generate histograms and derive geometric mean intensity readings for the binding of anti-CD27, -CD71 and CD86 antibodies to B cells. This data was imported into Excel and percentage inhibition values generated for each combination. The data was then imported into Graphpad Prism and box and whisker charts generated for each combination with the mean indicated by a '+'.
[0500] Figure 40 shows the inhibition of CD27 expression on B cells induced by VR4447/VR4126 BYbe and VR4447/VR4126/VR645 BYbe/Albumin. Across the five donors tested both showed consistently similar levels of inhibition of anti-IgM induced CD27. Figure 41 shows the inhibition of CD71 expression on B cells induced by VR4447/VR4126 BYbe and VR4447/VR4126/VR645 BYbe/Albumin. Across the five donors both showed consistently similar levels of inhibition of anti-IgM induced CD71. Figure 42 shows the inhibition of CD86 expression on B cells induced by VR4447/VR4126 BYbe and VR4447/VR4126/VR645 BYbe/Albumin. Across the five donors both showed consistently similar levels of inhibition of anti-IgM induced CD86.

GCN4(7P14P) sequences

ASGGGRMKQLEPKVEELLPKNYHLENEVARLKKLVGERHHHHHH SEQ ID NO: 1 wherein the amino acids in bold are optional and the amino acids in italics are the linking sequence

GCTAGCGGAGGCGGAAGAATGAAACAACCTTGAACCCAAGGTTGAAGAATTGCTTCCGAAAA
 TTATCACTTGGAAAATGAGGTTGCCAGATTAAAGAAATTAGTTGGCGAACGCCATCACCATC
 ACCATCAC SEQ ID NO: 2

52SR4 ds scFv sequence

DAVVTQESALTSSPGETVTTLTCRSSTGAVTTSNYASWVQEKPDHLFTGLIGGTNNRAPGV
 PARFSGSLIGDKAALTITGAQTEDEAIYFCVLWYSDHWVFGCGTKLTVLGGGGSGGGGS
 GGGSGGGGSDVQLQQSGPGLVAPSQSL SITCTVSGFLLTDYGVNWVRQSPGKCLEWLGV
 IWGDGITDYN SALKSRLSVTKDNSKSQVFLKMNSLQSGDSARYYCVTGLFDYWGGTTLT
 VSSAAAHHHHHHEQKLISEEDL— SEQ ID NO: 3

GATGCGGTGGTGACCCAGGAAAGCGCGCTGACCAGCAGCCCGGGCGAAACCGTGACCCTGAC
 CTGCCGCAGCAGACCCGGCGCGGTGACCACCAGCAACTATGCGAGCTGGGTGCAGGAAAAAC
 CGGATCATCTGTTTACCGGCCTGATTGGCGGCACCAACAACCGCGCGCCGGGCGTGCCGGCG
 CGCTTTAGCGGCAGCCTGATTGGCGATAAAGCGGCGCTGACCATTACCGGCGCGCAGACCGA
 AGATGAAGCGATTTATTTTTGCGTGCTGTGGTATAGCGACCATTGGGTGTTTGGCTGCGGCA
 CCAAACCTGACCGTGCTGGGTGGAGGCGGTGGCTCAGGCGGAGGTGGCTCAGGCGGTGGCGGG
 TCTGGCGGCGGCGGCAGCGATGTGCAGCTGCAGCAGAGCGGCCCGGGCCTGGTGGCGCCGAG
 CCAGAGCCTGAGCATTACCTGCACCGTGAGCGGCTTTCTCCTGACCGATTATGGCGTGAACT
 GGGTGCGCCAGAGCCCGGGCAAATGCCTGGAATGGCTGGGCGTGATTTGGGGCGATGGCATT
 ACCGATTATAACAGCGCGCTGAAAAGCCGCCTGAGCGTGACCAAAGATAACAGCAAAAGCCA
 GGTGTTTCTGAAAATGAACAGCCTGCAGAGCGGCGATAGCGCGCGCTATTATTGCGTGACCG
 GCCTGTTTGATTATTGGGGCCAGGGCACCAACCTGACCGTGAGCAGCGCGGCCGCCCATCAC
 CATCACCATCACGAACAGAACTGATTAGCGAAGAAGATCTGTAATAG
 SEQ ID NO: 4

SEQUENCE LISTING

[0501]

<110> UCB Pharma

<120> New Bispecific Format Suitable for Use in High-Through-Put Screening

<130> G0221_WO01

<160> 72

<170> Patent In version 3.5

<210> 1

<211> 44

<212> PRT

<213> Artificial Sequence

<220>

<223> GCN4(7P14P) sequence

<400> 1

Ala Ser Gly Gly Gly Arg Met Lys Gln Leu Glu Pro Lys Val Glu Glu
1 5 10 15

Leu Leu Pro Lys Asn Tyr His Leu Glu Asn Glu Val Ala Arg Leu Lys
20 25 30

10

Lys Leu Val Gly Glu Arg His His His His His His
35 40

15

<210>	2
<211>	132
<212>	DNA
<213>	Artificial Sequence

20 <220>
<223> DNA encoding SEQ ID NO: 1

<400> 2

25	gctagcggag gcggaagaat gaaacaactt gaacccaaggttgaagaatt gcttccgaaa	60
	aattatcact tggaaaatga gggtgccaga ttaaagaaat tagttggcga acgccatcac	120
	catcaccatc ac	132

30

<210> 3
<211> 262
<212> PRT
<213> Artificial Sequence

35 <220>
<223> 52SR4 ds scFv sequence

<400> 3

⁴⁰

Asp	Ala	Val	Val	Thr	Gln	Glu	Ser	Ala	Leu	Thr	Ser	Ser	Pro	Gly	Glu
1				5					10					15	

45

50

55

	Thr	Val	Thr	Leu	Thr	Cys	Arg	Ser	Ser	Thr	Gly	Ala	Val	Thr	Thr	Ser	
				20					25					30			
5	Asn	Tyr	Ala	Ser	Trp	Val	Gln	Glu	Lys	Pro	Asp	His	Leu	Phe	Thr	Gly	
			35					40					45				
10	Leu	Ile	Gly	Gly	Thr	Asn	Asn	Arg	Ala	Pro	Gly	Val	Pro	Ala	Arg	Phe	
		50					55					60					
15	Ser	Gly	Ser	Leu	Ile	Gly	Asp	Lys	Ala	Ala	Leu	Thr	Ile	Thr	Gly	Ala	
	65					70					75					80	
20	Gln	Thr	Glu	Asp	Glu	Ala	Ile	Tyr	Phe	Cys	Val	Leu	Trp	Tyr	Ser	Asp	
					85					90					95		
25	His	Trp	Val	Phe	Gly	Cys	Gly	Thr	Lys	Leu	Thr	Val	Leu	Gly	Gly	Gly	
				100					105					110			
30	Gly	Gly	Ser	Gly	Gly	Gly	Gly	Ser	Gly	Gly	Gly	Gly	Ser	Gly	Gly	Gly	
			115					120					125				
35	Gly	Ser	Asp	Val	Gln	Leu	Gln	Gln	Ser	Gly	Pro	Gly	Leu	Val	Ala	Pro	
		130					135					140					
40	Ser	Gln	Ser	Leu	Ser	Ile	Thr	Cys	Thr	Val	Ser	Gly	Phe	Leu	Leu	Thr	
	145					150					155					160	
45	Asp	Tyr	Gly	Val	Asn	Trp	Val	Arg	Gln	Ser	Pro	Gly	Lys	Cys	Leu	Glu	
					165					170					175		
50	Trp	Leu	Gly	Val	Ile	Trp	Gly	Asp	Gly	Ile	Thr	Asp	Tyr	Asn	Ser	Ala	
				180					185					190			
55	Leu	Lys	Ser	Arg	Leu	Ser	Val	Thr	Lys	Asp	Asn	Ser	Lys	Ser	Gln	Val	
			195					200					205				
60	Phe	Leu	Lys	Met	Asn	Ser	Leu	Gln	Ser	Gly	Asp	Ser	Ala	Arg	Tyr	Tyr	
		210					215					220					
65	Cys	Val	Thr	Gly	Leu	Phe	Asp	Tyr	Trp	Gly	Gln	Gly	Thr	Thr	Leu	Thr	
	225					230					235					240	
70	Val	Ser	Ser	Ala	Ala	Ala	His	His	His	His	His	His	Glu	Gln	Lys	Leu	
					245					250					255		
75	Ile	Ser	Glu	Glu	Asp	Leu											

260

5 <210> 4
 <211> 792
 <212> DNA
 <213> Artificial Sequence

10 <220>
 <223> DNA encoding SEQ ID NO: 3

 <400> 4

15	gatgcggtgg	tgaccagga	aagcgcgctg	accagcagcc	cgggcgaaac	cgtgaccctg	60
	acctgccgca	gcagcaccgg	cgcggtgacc	accagcaact	atgcgagctg	ggtgcaggaa	120
	aaaccggatc	atctgtttac	cggcctgatt	ggcggcacca	acaaccgcgc	gccgggcgctg	180
20	ccggcgcgct	ttagcggcag	cctgattggc	gataaagcgg	cgtgaccat	taccggcgcg	240
	cagaccgaag	atgaagcgat	ttatTTTTgc	gtgctgtggt	atagcgacca	ttgggtgttt	300
	ggctgcggca	ccaaactgac	cgtgctgggt	ggaggcggtg	gctcaggcgg	aggtggctca	360
25	ggcgggtggcg	ggtctggcgg	cggcggcagc	gatgtgcagc	tgacgcagag	cggcccgggc	420
	ctggtggcgc	cgagccagag	cctgagcatt	acctgcaccg	tgagcggctt	tctcctgacc	480
30	gattatggcg	tgaactgggt	gcgccagagc	ccgggcaaatt	gcctggaatg	gctgggcgctg	540
	atttggggcg	atggcattac	cgattataac	agcgcgctga	aaagccgcct	gagcgtgacc	600
	aaagataaca	gcaaaagcca	ggtgtttctg	aaaatgaaca	gcctgcagag	cggcgatagc	660
35	gcgcgctatt	attgcgtgac	cggcctgttt	gattattggg	gccagggcac	caccctgacc	720
	gtgagcagcg	cggccgcccc	tcaccatcac	catcacgaac	agaaactgat	tagcgaagaa	780
40	gatctgtaat	ag					792

45

<210>	5
<211>	8
<212>	PRT
<213>	Artificial Sequence

<220>
<223> FLEXIBLE PEPTIDE LINKER

50 <400> 5

Asp Lys Thr His Thr Cys Ala Ala
1 5

55 <210> 6
 <211> 11
 <212> PRT
 <213> Artificial Sequence

<220>

<223> FLEXIBLE PEPTIDE LINKER

<400> 6

5

Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala
1				5					10	

<210> 7

10

<211> 18

<212> PRT

<213> Artificial Sequence

<220>

15

<223> FLEXIBLE PEPTIDE LINKER

<400> 7

20

Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala	Thr	Cys	Pro	Pro	Cys
1				5				10					15		

Pro Ala

25

<210> 8

<211> 25

<212> PRT

<213> Artificial Sequence

30

<220>

<223> FLEXIBLE PEPTIDE LINKER

<400> 8

35

Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala	Thr	Cys	Pro	Pro	Cys
1				5				10					15		

Pro	Ala	Thr	Cys	Pro	Pro	Cys	Pro	Ala
			20				25	

40

<210> 9

<211> 30

<212> PRT

45

<213> Artificial Sequence

<220>

<223> FLEXIBLE PEPTIDE LINKER

50

<400> 9

Asp	Lys	Thr	His	Thr	Cys	Pro	Pro	Cys	Pro	Ala	Gly	Lys	Pro	Thr	Leu
1				5				10					15		

55

Tyr	Asn	Ser	Leu	Val	Met	Ser	Asp	Thr	Ala	Gly	Thr	Cys	Tyr
			20					25				30	

<210> 10
 <211> 31
 <212> PRT
 <213> Artificial Sequence
 5
 <220>
 <223> FLEXIBLE PEPTIDE LINKER
 <400> 10
 10
 Asp Lys Thr His Thr Cys Pro Pro Cys Pro Ala Gly Lys Pro Thr His
 1 5 10 15
 Val Asn Val Ser Val Val Met Ala Glu Val Asp Gly Thr Cys Tyr
 20 25 30
 <210> 11
 <211> 15
 <212> PRT
 <213> Artificial Sequence
 <220>
 <223> FLEXIBLE PEPTIDE LINKER
 25
 <400> 11
 Asp Lys Thr His Thr Cys Cys Val Glu Cys Pro Pro Cys Pro Ala
 1 5 10 15
 <210> 12
 <211> 26
 <212> PRT
 <213> Artificial Sequence
 35
 <220>
 <223> FLEXIBLE PEPTIDE LINKER
 <400> 12
 40
 Asp Lys Thr His Thr Cys Pro Arg Cys Pro Glu Pro Lys Ser Cys Asp
 1 5 10 15
 Thr Pro Pro Pro Cys Pro Arg Cys Pro Ala
 20 25
 <210> 13
 <211> 11
 <212> PRT
 <213> Artificial Sequence
 <220>
 <223> FLEXIBLE PEPTIDE LINKER
 55
 <400> 13

		Asp	Lys	Thr	His	Thr	Cys	Pro	Ser	Cys	Pro	Ala
		1				5					10	

5	<210> 14	
	<211> 7	
	<212> PRT	
	<213> Artificial Sequence	

10	<220>	
	<223> Flexible linker	

	<400> 14	
--	----------	--

15		Ser	Gly	Gly	Gly	Gly	Ser	Glu
		1				5		

	<210> 15	
	<211> 6	
	<212> PRT	
20	<213> Artificial Sequence	

	<220>	
	<223> Flexible linker	

25	<400> 15	
----	----------	--

		Asp	Lys	Thr	His	Thr	Ser
		1				5	

30	<210> 16	
	<211> 6	
	<212> PRT	
	<213> Artificial Sequence	

35	<220>	
	<223> Flexible linker	

	<400> 16	
--	----------	--

40		Ser	Gly	Gly	Gly	Gly	Ser
		1				5	

	<210> 17	
	<211> 11	
45	<212> PRT	
	<213> Artificial Sequence	

	<220>	
	<223> Flexible linker	

50	<400> 17	
----	----------	--

		Ser	Gly	Gly	Gly	Gly	Ser	Gly	Gly	Gly	Ser
		1			5					10	

55	<210> 18	
	<211> 16	
	<212> PRT	

<213> Artificial Sequence

<220>

<223> Flexible linker

5

<400> 18

10

Ser Gly Gly Gly Gly Ser Gly Gly Gly Gly Ser Gly Gly Gly Gly Ser
1 5 10 15

<210> 19

<211> 21

<212> PRT

<213> Artificial Sequence

15

<220>

<223> Flexible linker

<400> 19

20

Ser Gly Gly Gly Gly Ser Gly Gly Gly Gly Ser Gly Gly Gly Gly Ser
1 5 10 15

25

Gly Gly Gly Gly Ser
20

<210> 20

<211> 26

<212> PRT

<213> Artificial Sequence

30

<220>

<223> Flexible linker

35

<400> 20

Ser Gly Gly Gly Gly Ser Gly Gly Gly Gly Ser Gly Gly Gly Gly Ser
1 5 10 15

40

Gly Gly Gly Gly Ser Gly Gly Gly Gly Ser
20 25

45

<210> 21

<211> 11

<212> PRT

<213> Artificial Sequence

50

<220>

<223> Flexible linker

<400> 21

55

Ala Ala Ala Gly Ser Gly Gly Ala Ser Ala Ser
1 5 10

<210> 22

<211> 16
 <212> PRT
 <213> Artificial Sequence
 5 <220>
 <223> Flexible linker
 <220>
 <221> misc_feature
 10 <222> (7)..(7)
 <223> Xaa can be any naturally occurring amino acid
 <400> 22
 15 Ala Ala Ala Gly Ser Gly Xaa Gly Gly Gly Ser Gly Ala Ser Ala Ser
 1 5 10 15
 <210> 23
 <211> 21
 20 <212> PRT
 <213> Artificial Sequence
 <220>
 <223> Flexible linker
 25 <220>
 <221> misc_feature
 <222> (7)..(7)
 <223> Xaa can be any naturally occurring amino acid
 30 <220>
 <221> misc_feature
 <222> (12)..(12)
 <223> Xaa can be any naturally occurring amino acid
 35 <400> 23
 Ala Ala Ala Gly Ser Gly Xaa Gly Gly Gly Ser Xaa Gly Gly Gly Ser
 1 5 10 15
 40 Gly Ala Ser Ala Ser
 20
 45 <210> 24
 <211> 26
 <212> PRT
 <213> Artificial Sequence
 50 <220>
 <223> Flexible linker
 <220>
 <221> misc_feature
 55 <222> (7)..(7)
 <223> Xaa can be any naturally occurring amino acid
 <220>

<221> misc_feature
 <222> (12)..(12)
 <223> Xaa can be any naturally occurring amino acid
 5 <220>
 <221> misc_feature
 <222> (17)..(17)
 <223> Xaa can be any naturally occurring amino acid
 10 <400> 24

 Ala Ala Ala Gly Ser Gly Xaa Gly Gly Gly Ser Xaa Gly Gly Gly Ser
 1 5 10 15
 15
 Xaa Gly Gly Gly Ser Gly Ala Ser Ala Ser
 20 25

 <210> 25
 20 <211> 31
 <212> PRT
 <213> Artificial Sequence

 <220>
 25 <223> Flexible linker

 <220>
 <221> misc_feature
 <222> (7)..(7)
 30 <223> Xaa can be any naturally occurring amino acid

 <220>
 <221> misc_feature
 <222> (12)..(12)
 35 <223> Xaa can be any naturally occurring amino acid

 <220>
 <221> misc_feature
 <222> (17)..(17)
 40 <223> Xaa can be any naturally occurring amino acid

 <220>
 <221> misc_feature
 <222> (22)..(22)
 45 <223> Xaa can be any naturally occurring amino acid

 <400> 25

 Ala Ala Ala Gly Ser Gly Xaa Gly Gly Gly Ser Xaa Gly Gly Gly Ser
 1 5 10 15
 50
 Xaa Gly Gly Gly Ser Xaa Gly Gly Gly Ser Gly Ala Ser Ala Ser
 20 25 30
 55
 <210> 26
 <211> 13
 <212> PRT

<213> Artificial Sequence

<220>

<223> Flexible linker

<220>

<221> misc_feature

<222> (7)..(7)

<223> Xaa can be any naturally occurring amino acid

<400> 26

Ala Ala Ala Gly Ser Gly Xaa Ser Gly Ala Ser Ala Ser
1 5 10

<210> 27

<211> 28

<212> PRT

<213> Artificial Sequence

<220>

<223> Flexible linker

<400> 27

Pro Gly Gly Asn Arg Gly Thr Thr Thr Thr Arg Arg Pro Ala Thr Thr
1 5 10 15

Thr Gly Ser Ser Pro Gly Pro Thr Gln Ser His Tyr
20 25

<210> 28

<211> 11

<212> PRT

<213> Artificial Sequence

<220>

<223> Flexible linker

<400> 28

Ala Thr Thr Thr Gly Ser Ser Pro Gly Pro Thr
1 5 10

<210> 29

<211> 6

<212> PRT

<213> Artificial Sequence

<220>

<223> Flexible linker

<400> 29

Ala Thr Thr Thr Gly Ser
1 5

<210> 30
 <211> 13
 <212> PRT
 <213> Artificial Sequence
 5
 <220>
 <223> Flexible linker
 <400> 30
 10
 Ala Ala Ala Ala Ala Ala Ala Ala Ala Ala Ala Ala Ala Ala
 1 5 10
 <210> 31
 <211> 21
 <212> PRT
 <213> Artificial Sequence
 15
 <220>
 <223> Flexible linker
 20
 <400> 31
 25
 Glu Pro Ser Gly Pro Ile Ser Thr Ile Asn Ser Pro Pro Ser Lys Glu
 1 5 10 15
 Ser His Lys Ser Pro
 20
 30
 <210> 32
 <211> 15
 <212> PRT
 <213> Artificial Sequence
 35
 <220>
 <223> Flexible linker
 <400> 32
 40
 Gly Thr Val Ala Ala Pro Ser Val Phe Ile Phe Pro Pro Ser Asp
 1 5 10 15
 <210> 33
 <211> 15
 <212> PRT
 <213> Artificial Sequence
 45
 <220>
 <223> Flexible linker
 50
 <400> 33
 55
 Gly Gly Gly Gly Ile Ala Pro Ser Met Val Gly Gly Gly Gly Ser
 1 5 10 15
 <210> 34
 <211> 15

<212> PRT
 <213> Artificial Sequence

 <220>
 5 <223> Flexible linker

 <400> 34

 10 Gly Gly Gly Gly Lys Val Glu Gly Ala Gly Gly Gly Gly Gly Ser
 1 5 10 15

 <210> 35
 <211> 15
 <212> PRT
 15 <213> Artificial Sequence

 <220>
 <223> Flexible linker

 20 <400> 35

 Gly Gly Gly Gly Ser Met Lys Ser His Asp Gly Gly Gly Gly Ser
 1 5 10 15

 25 <210> 36
 <211> 15
 <212> PRT
 <213> Artificial Sequence

 30 <220>
 <223> Flexible linker

 <400> 36

 35 Gly Gly Gly Gly Asn Leu Ile Thr Ile Val Gly Gly Gly Gly Ser
 1 5 10 15

 <210> 37
 <211> 15
 40 <212> PRT
 <213> Artificial Sequence

 <220>
 <223> Flexible linker
 45
 <400> 37

 Gly Gly Gly Gly Val Val Pro Ser Leu Pro Gly Gly Gly Gly Ser
 1 5 10 15
 50
 <210> 38
 <211> 12
 <212> PRT
 <213> Artificial Sequence
 55
 <220>
 <223> Flexible linker

<400> 38

5

Gly Gly Glu Lys Ser Ile Pro Gly Gly Gly Gly Ser
 1 5 10

<210> 39

<211> 18

<212> PRT

<213> Artificial Sequence

10

<220>

<223> Flexible linker

<400> 39

15

Arg Pro Leu Ser Tyr Arg Pro Pro Phe Pro Phe Gly Phe Pro Ser Val
 1 5 10 15

20

Arg Pro

<210> 40

<211> 18

<212> PRT

<213> Artificial Sequence

25

<220>

<223> Flexible linker

30

<400> 40

Tyr Pro Arg Ser Ile Tyr Ile Arg Arg Arg His Pro Ser Pro Ser Leu
 1 5 10 15

35

Thr Thr

<210> 41

<211> 18

<212> PRT

<213> Artificial Sequence

40

<220>

<223> Flexible linker

45

<400> 41

Thr Pro Ser His Leu Ser His Ile Leu Pro Ser Phe Gly Leu Pro Thr
 1 5 10 15

50

Phe Asn

<210> 42

<211> 18

<212> PRT

<213> Artificial Sequence

55

<220>

<223> Flexible linker

<400> 42

5

Arg	Pro	Val	Ser	Pro	Phe	Thr	Phe	Pro	Arg	Leu	Ser	Asn	Ser	Trp	Leu
1				5					10					15	

10

Pro Ala

<210> 43

<211> 18

<212> PRT

15

<213> Artificial Sequence

<220>

<223> Flexible linker

20

<400> 43

Ser	Pro	Ala	Ala	His	Phe	Pro	Arg	Ser	Ile	Pro	Arg	Pro	Gly	Pro	Ile
1				5					10					15	

25

Arg Thr

<210> 44

<211> 18

<212> PRT

30

<213> Artificial Sequence

<220>

<223> Flexible linker

35

<400> 44

Ala	Pro	Gly	Pro	Ser	Ala	Pro	Ser	His	Arg	Ser	Leu	Pro	Ser	Arg	Ala
1				5					10					15	

40

Phe Gly

<210> 45

<211> 18

<212> PRT

45

<213> Artificial Sequence

<220>

<223> Flexible linker

50

<400> 45

Pro	Arg	Asn	Ser	Ile	His	Phe	Leu	His	Pro	Leu	Leu	Val	Ala	Pro	Leu
1				5					10					15	

55

Gly Ala

<210> 46
 <211> 18
 <212> PRT
 <213> Artificial Sequence
 5
 <220>
 <223> Flexible linker
 10
 <400> 46
 Met Pro Ser Leu Ser Gly Val Leu Gln Val Arg Tyr Leu Ser Pro Pro
 1 5 10 15
 15 Asp Leu
 <210> 47
 <211> 18
 <212> PRT
 20 <213> Artificial Sequence
 <220>
 <223> Flexible linker
 25 <400> 47
 Ser Pro Gln Tyr Pro Ser Pro Leu Thr Leu Thr Leu Pro Pro His Pro
 1 5 10 15
 30 Ser Leu
 <210> 48
 <211> 18
 35 <212> PRT
 <213> Artificial Sequence
 <220>
 <223> Flexible linker
 40 <400> 48
 Asn Pro Ser Leu Asn Pro Pro Ser Tyr Leu His Arg Ala Pro Ser Arg
 1 5 10 15
 45 Ile Ser
 <210> 49
 50 <211> 17
 <212> PRT
 <213> Artificial Sequence
 <220>
 55 <223> Flexible linker
 <400> 49

Leu Pro Trp Arg Thr Ser Leu Leu Pro Ser Leu Pro Leu Arg Arg Arg
 1 5 10 15

5 Pro

10 <210> 50
 <211> 18
 <212> PRT
 <213> Artificial Sequence

15 <220>
 <223> Flexible linker

<400> 50

20 Pro Pro Leu Phe Ala Lys Gly Pro Val Gly Leu Leu Ser Arg Ser Phe
 1 5 10 15

Pro Pro

25 <210> 51
 <211> 18
 <212> PRT
 <213> Artificial Sequence

30 <220>
 <223> Flexible linker

<400> 51

35 Val Pro Pro Ala Pro Val Val Ser Leu Arg Ser Ala His Ala Arg Pro
 1 5 10 15

Pro Tyr

40 <210> 52
 <211> 17
 <212> PRT
 <213> Artificial Sequence

45 <220>
 <223> Flexible linker

<400> 52

50 Leu Arg Pro Thr Pro Pro Arg Val Arg Ser Tyr Thr Cys Cys Pro Thr
 1 5 10 15

55 Pro

<210> 53
 <211> 18

<212> PRT
 <213> Artificial Sequence

<220>
 <223> Flexible linker

<400> 53

10 Pro Asn Val Ala His Val Leu Pro Leu Leu Thr Val Pro Trp Asp Asn
 1 5 10 15

Leu Arg

15 <210> 54
 <211> 18
 <212> PRT
 <213> Artificial Sequence

20 <220>
 <223> Flexible linker

<400> 54

25 Cys Asn Pro Leu Leu Pro Leu Cys Ala Arg Ser Pro Ala Val Arg Thr
 1 5 10 15

Phe Pro

30 <210> 55
 <211> 11
 <212> PRT
 <213> Artificial Sequence

35 <220>
 <223> Linker

<400> 55

40 Asp Leu Cys Leu Arg Asp Trp Gly Cys Leu Trp
 1 5 10

45 <210> 56
 <211> 11
 <212> PRT
 <213> Artificial Sequence

50 <220>
 <223> Linker

<400> 56

55 Asp Ile Cys Leu Pro Arg Trp Gly Cys Leu Trp
 1 5 10

<210> 57
 <211> 15

<212> PRT
 <213> Artificial Sequence

 <220>
 5 <223> Linker

 <400> 57

 10 Met Glu Asp Ile Cys Leu Pro Arg Trp Gly Cys Leu Trp Gly Asp
 1 5 10 15

 <210> 58
 <211> 20
 <212> PRT
 15 <213> Artificial Sequence

 <220>
 <223> Linker

 20 <400> 58

 Gln Arg Leu Met Glu Asp Ile Cys Leu Pro Arg Trp Gly Cys Leu Trp
 1 5 10 15

 25 Glu Asp Asp Glu
 20

 <210> 59
 30 <211> 20
 <212> PRT
 <213> Artificial Sequence

 <220>
 35 <223> Linker

 <400> 59

 40 Gln Gly Leu Ile Gly Asp Ile Cys Leu Pro Arg Trp Gly Cys Leu Trp
 1 5 10 15

 Gly Arg Ser Val
 20

 45 <210> 60
 <211> 21
 <212> PRT
 <213> Artificial Sequence
 50
 <220>
 <223> Linker

 <400> 60
 55

Gln Gly Leu Ile Gly Asp Ile Cys Leu Pro Arg Trp Gly Cys Leu Trp
 1 5 10 15

5 Gly Arg Ser Val Lys
 20

<210> 61
 <211> 15
 10 <212> PRT
 <213> Artificial Sequence

<220>
 <223> Linker

15 <400> 61

Glu Asp Ile Cys Leu Pro Arg Trp Gly Cys Leu Trp Glu Asp Asp
 1 5 10 15

20 <210> 62
 <211> 18
 <212> PRT
 <213> Artificial Sequence

25 <220>
 <223> Linker

<400> 62

30 Arg Leu Met Glu Asp Ile Cys Leu Pro Arg Trp Gly Cys Leu Trp Glu
 1 5 10 15

35 Asp Asp

<210> 63
 <211> 16
 <212> PRT
 40 <213> Artificial Sequence

<220>
 <223> Linker

45 <400> 63

Met Glu Asp Ile Cys Leu Pro Arg Trp Gly Cys Leu Trp Glu Asp Asp
 1 5 10 15

50 <210> 64
 <211> 15
 <212> PRT
 <213> Artificial Sequence

55 <220>
 <223> Linker

<400> 64

		Met	Glu	Asp	Ile	Cys	Leu	Pro	Arg	Trp	Gly	Cys	Leu	Trp	Glu	Asp	
		1				5					10					15	
5		<210> 65															
		<211> 18															
		<212> PRT															
		<213> Artificial Sequence															
		<220>															
10		<223> Linker															
		<400> 65															
15		Arg	Leu	Met	Glu	Asp	Ile	Cys	Leu	Ala	Arg	Trp	Gly	Cys	Leu	Trp	Glu
		1				5					10					15	
		Asp Asp															
20		<210> 66															
		<211> 20															
		<212> PRT															
		<213> Artificial Sequence															
		<220>															
25		<223> Linker															
		<400> 66															
30		Glu	Val	Arg	Ser	Phe	Cys	Thr	Arg	Trp	Pro	Ala	Glu	Lys	Ser	Cys	Lys
		1				5					10					15	
		Pro Leu Arg Gly															
35						20											
		<210> 67															
		<211> 20															
		<212> PRT															
40		<213> Artificial Sequence															
		<220>															
		<223> Linker															
45		<400> 67															
		Arg	Ala	Pro	Glu	Ser	Phe	Val	Cys	Tyr	Trp	Glu	Thr	Ile	Cys	Phe	Glu
		1				5					10					15	
50		Arg Ser Glu Gln															
						20											
		<210> 68															
55		<211> 11															
		<212> PRT															
		<213> Artificial Sequence															

<220>
 <223> Linker

 <400> 68
 5

 Glu Met Cys Tyr Phe Pro Gly Ile Cys Trp Met

 10 1 5 10

 <210> 69
 <211> 11
 <212> PRT
 15 <213> Artificial Sequence

 <220>
 <223> peptide sequence

 20 <400> 69

 Gly Ala Pro Ala Pro Ala Ala Pro Ala Pro Ala
 1 5 10

 25 <210> 70
 <211> 4
 <212> PRT
 <213> Artificial Sequence

 30 <220>
 <223> peptide sequence

 <400> 70

 35 Pro Pro Pro Pro
 1

 <210> 71
 <211> 6
 40 <212> PRT
 <213> Artificial Sequence

 <220>
 <223> Linker
 45
 <400> 71

 Ala Ser Gly Gly Gly Gly
 1 5
 50
 <210> 72
 <211> 8
 <212> PRT
 <213> Artificial Sequence
 55
 <220>
 <223> Linker

<400> 72

5

Ala Ser Gly Gly Gly Gly Ser Gly
1 5

REFERENCES CITED IN THE DESCRIPTION

This list of references cited by the applicant is for the reader's convenience only. It does not form part of the European patent document. Even though great care has been taken in compiling the references, errors or omissions cannot be excluded and the EPO disclaims all liability in this regard.

Patent documents cited in the description

- WO 2014001326 A [0021]
- US 5731168 A [0025]
- US 7183076 B [0025]
- WO 9311162 A [0027]
- WO 05003169 A [0171]
- WO 05003170 A [0171]
- WO 05003171 A [0171]
- WO 9222853 A [0171]
- WO 05113605 A [0171]
- WO 2009040562 A [0171] [0251]
- WO 2010035012 A [0171] [0251]
- WO 9820734 A [0324]
- EP 1570267 B1 [0433]

Non-patent literature cited in the description

- D. HOLMES. *Nature Rev Drug Disc*, November 2011, vol. 10, 798 [0005]
- RIDGWAY et al. *Protein Eng.*, 1996, vol. 9, 617-621 [0025]
- MERCHANT et al. *Nat. Biotechnol.*, 1998, vol. 16 (7), 677-81 [0025]
- HOLLIGER ; HUDSON. *Nature Biotech.*, 2005, vol. 23 (9), 1126-1136 [0171]
- ADAIR ; LAWSON. *Drug Design Reviews - Online*, 2005, vol. 2 (3), 209-217 [0171]
- VERMA et al. *Journal of Immunological Methods*, 1998, vol. 216, 165-181 [0171]
- YANG et al. *J. Mol. Biol.*, 1995, vol. 254, 392-403 [0201]
- MARKS et al. *Bio/Technology*, 1992, vol. 10, 779-783 [0201]
- LOW et al. *J. Mol. Biol.*, 1996, vol. 250, 359-368 [0201]
- PATTEN et al. *Curr. Opin. Biotechnol.*, 1997, vol. 8, 724-733 [0201]
- THOMPSON et al. *J. Mol. Biol.*, 1996, vol. 256, 77-88 [0201]
- CRAMERI et al. *Nature*, 1998, vol. 391, 288-291 [0201]
- Computational Molecular Biology. Oxford University Press, 1988 [0226]
- Biocomputing. Informatics and Genome Projects. Academic Press, 1993 [0226]
- Computer Analysis of Sequence Data. Humana Press, 1994 [0226]
- HEINJE, G. *Sequence Analysis in Molecular Biology*. Academic Press, 1987 [0226]
- *Sequence Analysis Primer*. M Stockton Press, 1991 [0226]
- ALTSCHUL, S.F. et al. *J. Mol. Biol.*, 1990, vol. 215, 403-410 [0226]
- GISH, W. ; STATES, D.J. *Nature Genet.*, 1993, vol. 3, 266-272 [0226]
- MADDEN, T.L. et al. *Meth. Enzymol.*, 1996, vol. 266, 131-141 [0226]
- ALTSCHUL, S.F. et al. *Nucleic Acids Res.*, 1997, vol. 25, 3389-3402 [0226]
- ZHANG, J. ; MADDEN, T.L. *Genome Res.*, 1997, vol. 7, 649-656 [0226]
- ANGAL et al. *Molecular Immunology*, 1993, vol. 30, 105-108 [0275]
- HARRIS, R.J. *Journal of Chromatography*, 1995, vol. 705, 129-134 [0276]
- Remington's Pharmaceutical Sciences. Mack Publishing Company, 1991 [0313]
- Current Protocols in Molecular Biology. Wiley Interscience, 1999 [0347]
- Maniatis Manual. Cold Spring Harbor Publishing [0347]
- HANES J ; JERMUTUS L ; WEBER-BORNHAUSER S ; BOSSHARD HR ; PLÜCKTHUN A. *Proc. Natl. Acad. Sci. U.S.A.*, 1998, vol. 95, 14130-14135 [0363]
- ZHAND C ; SPINELLI S ; LUGINBUHL B ; AMSTUTZ P ; CAMBILLAU C ; PLUCKTHUN A. *J. Biol. Chem.*, 2004, vol. 279, 18870-18877 [0363]
- BERGER C ; WEBER-BORNHAUSER S ; EGGENBERGER Y ; HANES J ; PLUCKTHUN A ; BOSSHARD H. R. *F.E.B.S. Letters*, 1999, vol. 450, 149-153 [0363]
- CLARGO et al. *Mabs*, 01 January 2014, vol. 6 (1), 143-159 [0433]

PATENTKRAV

1. Bispecifikt proteinkompleks med formlen A-X:Y-B hvor:

A-X er et første fusionsprotein;

5 Y-B er et andet fusionsprotein;

X:Y er et heterodimert bånd;

: er en bindingsinteraktion mellem X og Y;

A er en første proteinbestanddel af det bispecifikke proteinkompleks valgt blandt et Fab- eller Fab' fragment;

10 B er en anden proteinbestanddel af det bispecifikke proteinkompleks valgt blandt et Fab- eller Fab'-fragment;

hvor X uafhængigt er valgt blandt et scFv eller et peptid, forudsat, når X er et peptid, Y er et scFv, og når X er scFv, Y er et peptid;

hvor X eller Y er et scFv specifikt for peptidet GCN4 repræsenteret ved SEQ
15 ID NO:1 eller aminosyrerne 1 til 38 ifølge SEQ ID NO: 1; hvor scFv'et betegnes 52SR4 og er repræsenteret ved SEQ ID NO: 3 eller aminosyrerne 1 til 243 ifølge SEQ ID NO: 3, og hvor X eller Y er et peptid GCN4 repræsenteret ved SEQ ID NO: 1 eller aminosyrerne 1 til 38 ifølge SEQ ID NO: 1.

2. Bispecifikt proteinkompleks ifølge krav 1, hvor A er et Fab-fragment;
20 og/eller B er et Fab-fragment.

3. Bispecifikt proteinkompleks ifølge et hvilket som helst af kravene 1 eller 2, hvor X og/eller Y er fusioneret, eventuelt via en linker, til C-terminalen af tungkæden i Fab- eller Fab'-fragmentet.

4. Bispecifikt proteinkompleks ifølge et hvilket som helst af kravene 1 til 3,
25 hvor A og/eller B er specifikt for et antigen valgt fra gruppen omfattende: celleoverfladereceptorer, såsom T-celle- eller B-cellessignaleringsreceptorer, co-stimulerende molekyler, checkpoint-hæmmer, naturlige dræbercellereceptorer, Immunglobulinreceptorer, TNFR-familiereceptorer, B7-familiereceptorer, adhæsionsmolekyler, integriner, cytokin-/chemokinreceptorer, GPCR,


vækstfaktorreceptorer, kinasereceptorer, vævsspecifikke antigener, cancerantigener, patogengenkendelsesreceptorer, komplementreceptorer, hormonreceptorer eller opløselige molekyler, såsom cytokiner, chemokiner, leukotriener, vækstfaktorer, hormoner, enzymer og ionkanaler.

5 **5.** Sammensætning, der omfatter ét eller flere bispecifikke proteinkomplekser defineret i et hvilket som helst af kravene 1 til 4.

6. Bispecifikt proteinkompleks ifølge et hvilket som helst af kravene 1 til 4 eller sammensætning ifølge krav 5 til anvendelse i terapi.

DRAWINGS

Figure 1






Figure 2

	A₁X	A₂X	A₃X	A₄X
B₁Y	A ₁ X B ₁ Y	A ₂ X B ₁ Y	A ₃ X B ₁ Y	A ₄ X B ₁ Y
B₂Y	A ₁ X B ₂ Y	A ₂ X B ₂ Y	A ₃ X B ₂ Y	A ₄ X B ₂ Y
B₃Y	A ₁ X B ₃ Y	A ₂ X B ₃ Y	A ₃ X B ₃ Y	A ₄ X B ₃ Y
B₄Y	A ₁ X B ₄ Y	A ₂ X B ₄ Y	A ₃ X B ₄ Y	A ₄ X B ₄ Y

Figure 3




Figure 4





Figure 5

Figure 6

Surface plasmon resonance analysis of Antigen5Fab- scFv (52RS4) binding to GNC4 peptide

Sample	Biotin-GCN4 Peptide			
	ka (1/Ms)	kd (1/s)	Affinity (pM)	Immuno-reactivity %
Antigen5Fab52RS4scFv	5.10E+05	2.70E-04	516	105

% Immunoreactivity = $R_{\text{max}}/\text{Estimated } R_{\text{max}} \times 100$






Figure 7




Figure 8




Figure 9




Figure 10




Figure 11




Figure 12




Figure 13




Figure 14




Figure 15




Figure 16




Figure 17




Figure 18




Figure 19




Figure 20

Figure 21

Fab-X specificity antigen 1 -

SP_Y...	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
1	4.37	33.57	10.45	13.99	3.81	4.26	8.34	21.07	4.10	13.98	6.07	16.03	7.39	9.00	23.01
2	27.22	25.26	52.32	63.81	30.98	19.94	28.90	29.48	20.58	30.61	9.05	4.18	14.73	36.56	33.04
3	17.62	69.66	3.88	16.35	7.81	-5.45	6.94	12.43	6.90	7.41	15.55	18.23	8.31	7.13	15.98
4	13.81	70.47	-9.21	11.98	2.91	-0.27	7.92	9.85	0.93	4.02	8.95	40.06	10.33	6.02	15.96
5	5.05	36.09	-1.78	-0.25	10.55	-14.55	11.68	6.91	-0.32	-2.82	9.14	12.00	3.48	7.99	9.03
6	---	21.62	7.24	15.78	16.88	2.89	2.50	14.67	9.48	4.44	11.14	20.93	12.37	---	---
7	10.54	33.72	2.40	16.16	14.92	1.39	2.23	14.14	5.84	10.62	8.02	17.58	9.49	10.65	10.80
8	6.03	26.05	-7.89	3.46	7.85	-18.33	0.26	-5.98	-2.08	0.09	7.83	5.13	2.60	9.59	9.47
9	13.56	32.14	-4.74	2.87	9.63	-7.13	5.17	10.72	0.79	1.69	9.22	13.53	2.49	15.70	10.52
10	15.35	54.98	13.21	17.70	-2.07	6.41	21.65	23.25	9.94	10.33	16.30	41.50	13.23	6.52	16.72
11	---	-0.28	-50.76	-5.43	-0.23	-118.10	-4.78	-8.59	0.08	-4.95	6.09	2.91	-9.06	---	---
12	5.85	8.47	19.70	42.01	20.13	12.42	15.14	26.89	20.89	24.69	18.86	36.32	18.07	20.28	38.97
13	---	31.56	22.19	28.34	15.50	13.12	14.23	27.59	22.46	14.25	16.00	26.84	16.21	---	---
14	9.13	40.70	10.80	14.56	5.88	8.89	7.52	22.68	6.27	9.04	7.34	33.39	11.20	7.01	14.84
15	9.38	33.18	7.13	11.20	15.82	11.35	9.90	19.61	9.90	10.62	13.98	28.99	15.26	13.77	25.64
16	---	26.46	-7.59	19.38	15.65	-24.93	12.05	14.00	15.27	11.98	13.29	20.57	11.99	---	---
17	---	22.79	-4.70	19.07	17.98	-37.56	10.99	11.71	15.83	8.81	10.98	21.28	9.71	---	---
18	---	28.91	29.09	33.25	23.18	10.61	13.85	26.28	21.37	18.39	18.45	32.78	18.27	---	---
19	---	33.63	-15.93	2.00	15.54	-24.55	11.73	10.55	7.32	4.45	15.19	20.98	10.61	---	---
20	---	24.94	-13.88	20.88	10.53	-21.56	8.43	6.69	8.81	7.57	10.73	17.52	13.55	---	---
21	---	23.94	6.69	22.70	22.45	-3.18	13.33	21.28	15.88	10.81	15.92	25.96	17.22	---	---
22	13.93	20.55	-18.75	-4.96	6.90	-17.32	4.12	-3.15	-4.19	-6.48	11.71	2.48	-0.84	9.14	11.68
23	---	24.25	-9.45	22.83	24.63	-26.70	9.65	13.79	10.74	7.96	13.02	24.43	18.33	---	---

% Inhibition Syk

Figure 22

Fab-X specificity antigen 1 -15

SP_Y...	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
1	0.81	22.29	11.45	10.78	5.39	1.60	7.43	5.90	2.24	3.29	-2.71	10.80	0.99	6.90	-4.08
2	19.70	18.61	44.00	50.83	23.17	19.01	24.57	21.47	25.89	39.80	12.35	0.52	18.91	28.55	25.89
3	17.93	61.18	3.00	15.38	5.38	-0.12	4.05	3.77	6.84	12.78	11.01	16.71	2.97	9.17	4.69
4	13.34	55.83	-3.82	7.24	1.41	-1.39	6.75	0.89	3.12	3.12	4.16	40.49	1.71	5.13	1.10
5	0.34	25.16	0.51	-0.63	7.36	-1.41	11.49	4.03	5.59	2.99	5.76	11.67	2.84	5.90	-2.06
6	---	20.71	-0.85	6.04	0.40	3.49	-2.60	-0.87	3.41	0.22	1.34	12.34	3.36	---	---
7	10.95	27.95	1.38	15.66	11.27	2.40	2.48	4.45	6.13	16.07	3.18	12.77	3.92	5.85	0.37
8	7.73	24.86	4.02	5.23	3.56	0.49	2.68	-1.16	2.39	1.45	4.10	8.50	-0.64	3.01	0.97
9	7.37	28.27	-0.02	3.22	6.39	1.35	4.75	4.06	4.92	4.09	3.10	10.50	0.62	5.62	1.50
10	17.08	47.69	7.72	5.65	-5.55	1.40	21.89	6.63	5.90	5.11	7.61	35.77	1.50	4.32	-0.06
11	---	11.76	-10.93	1.38	-0.56	-24.94	-1.36	0.83	3.58	3.88	5.92	11.18	-4.57	---	---
12	6.05	-3.88	9.02	29.48	4.92	1.15	1.89	5.75	7.23	20.95	5.44	17.92	3.37	16.34	8.67
13	---	23.78	-0.99	7.30	0.43	-1.40	1.83	6.41	8.91	3.25	2.76	11.64	5.41	---	---
14	4.40	28.11	1.77	5.61	4.01	1.73	5.95	2.45	4.05	5.85	-0.60	20.60	5.45	4.11	0.05
15	4.96	25.60	0.87	3.43	5.41	1.77	1.89	1.89	3.72	1.96	2.29	14.57	5.26	2.36	0.11
16	---	22.39	-5.77	4.50	2.41	-6.86	4.89	2.52	5.19	3.07	3.38	9.47	1.17	---	---
17	---	19.53	-6.98	6.73	2.99	-10.27	2.35	2.13	6.02	2.08	3.10	10.86	2.03	---	---
18	---	19.90	2.43	12.38	2.08	-1.05	1.09	4.25	6.70	4.90	4.40	12.49	4.34	---	---
19	---	35.50	-12.22	-2.38	-1.61	-9.05	5.97	1.78	5.27	4.63	5.92	15.60	3.16	---	---
20	---	21.53	-9.67	2.95	-0.16	-6.71	-3.46	-4.81	4.13	2.31	6.52	10.96	4.25	---	---
21	---	19.03	-6.66	2.36	1.10	-6.49	-0.95	5.09	4.96	3.55	3.97	10.89	3.63	---	---
22	14.73	26.47	-5.00	0.58	3.97	-0.08	3.70	-1.25	1.72	0.99	11.40	7.99	2.80	8.39	4.58
23	---	19.72	-16.22	0.75	-1.16	-11.03	-3.57	3.30	3.57	1.23	4.90	8.22	4.20	---	---

% Inhibition

Fab-Y specificity antigen 1 -23

Figure 23

Fab-X specificity antigen 1-15

SP_Y...	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
1	-11.06	53.06	-1.56	6.02	-8.08	-12.41	3.94	-12.93	-8.22	-8.88	-11.49	-3.34	-17.29	-4.43	-6.28
2	40.41	37.15	48.29	48.35	39.50	30.66	40.36	25.61	42.23	42.64	32.59	37.11	26.58	35.70	38.11
3	13.19	70.81	2.37	15.53	8.68	7.93	12.32	3.78	10.56	14.18	8.76	-1.58	6.61	10.75	10.83
4	8.35	55.50	-6.07	17.80	7.97	7.14	8.49	6.54	6.98	9.87	12.96	41.04	5.99	9.28	7.89
5	-3.80	40.21	-0.85	13.54	2.05	-2.13	3.01	-6.86	0.91	2.02	0.76	4.48	-3.42	4.06	0.15
6	---	40.39	6.67	19.26	-4.86	12.95	3.52	2.59	6.73	5.65	5.28	4.77	-1.03	---	---
7	0.53	45.36	4.87	18.76	4.08	0.72	-1.43	5.07	4.69	11.49	6.48	5.72	4.78	0.73	-3.71
8	6.61	37.54	5.72	16.25	-1.45	2.38	0.64	0.17	1.10	2.08	6.50	7.58	2.18	-4.37	-4.53
9	0.51	42.34	0.95	9.41	-0.30	1.07	2.86	-2.56	0.51	-0.55	2.07	2.49	-2.84	4.54	0.16
10	11.13	52.32	3.75	12.67	-4.41	3.82	15.78	-3.15	5.70	4.26	1.77	22.73	-4.68	7.65	8.05
11	---	33.40	-3.95	18.04	7.48	6.00	3.13	4.24	5.66	13.11	9.85	11.24	2.04	---	---
12	0.69	29.50	-7.93	38.94	2.68	4.95	-0.13	-2.53	2.52	15.03	-0.82	3.63	-5.61	25.99	12.57
13	---	35.33	-4.80	10.35	-4.98	-7.45	-7.91	-6.47	3.78	-0.83	-4.68	-3.84	-9.54	---	---
14	1.62	44.02	0.91	16.05	1.80	-2.56	1.09	-4.49	3.94	5.30	-0.19	8.15	-5.85	5.21	0.66
15	2.25	36.02	-2.44	10.00	-0.22	0.37	-2.20	-2.11	2.41	1.84	3.23	2.78	2.51	-2.19	-4.62
16	---	33.56	-3.12	15.33	0.19	-5.46	-0.71	-5.67	1.83	1.35	-0.24	-1.99	-8.53	---	---
17	---	35.02	-4.92	16.87	1.78	-2.20	0.34	-2.12	3.07	1.17	0.98	0.80	-2.14	---	---
18	---	31.31	-1.77	19.71	-6.19	-4.39	-3.72	-3.91	2.70	3.44	-0.42	-1.65	-5.72	---	---
19	---	50.00	-2.93	9.92	2.01	9.54	10.19	2.60	11.44	9.84	7.96	9.17	3.37	---	---
20	---	38.56	-1.83	18.96	1.40	5.57	0.37	1.23	4.23	6.85	10.48	6.54	2.60	---	---
21	---	36.74	1.67	11.59	-5.60	6.23	0.59	6.00	6.99	5.11	3.12	3.73	3.44	---	---
22	-2.71	40.50	-2.45	12.91	-1.39	6.77	1.56	4.97	1.68	2.93	18.10	8.41	8.08	-2.35	-4.93
23	---	33.54	-1.68	5.45	-2.46	7.13	0.05	-3.14	2.99	3.34	6.02	4.97	-4.12	---	---

% Inhibition AKT

Fab-Y specificity antigen 1-23




Figure 24




Figure 25




Figure 26




Figure 27

% Inhibition (Mean $\pm$ SD) of CD71 on Anti-IgM Stimulated B-cells




Figure 28

% Inhibition (Mean $\pm$ SD) of Phosphorylated p38 in Anti-IgM Stimulated B-cells




Figure 29

Inhibition of phospho readouts by purified Antigen 2-specific Fab-Y (VR4447) combined with Antigen 3-specific Fab-X transients on IgM stimulated B-cells from UCB_Cone_172




Figure 30




Figure 31




Figure 32




Figure 33

Figure 34

**% Inhibition (+/- SD) of Phosphorylated PLC γ 2 in
Anti-IgM Stimulated B-cells
by Antigen 3 and Antigen 2 Specific Fab-Kd-Fab and BYbe.**




Figure 35

**% Inhibition (+/- SD) of Phosphorylated P38 in
Anti-IgM Stimulated B-cells,
by Antigen 3 and Antigen 2 Specific Fab-Kd-Fab and BYbe.**




Figure 36

**% Inhibition (+/- SD) of Phosphorylated Akt in
Anti-IgM Stimulated B-cells,
by Antigen 3 and Antigen 2 Specific Fab-Kd-Fab and BYbe.**




Figure 37

**% Inhibition of Anti-IgM Induced CD71 on B-cells,
by Antigen 3 and Antigen 2 Specific Fab-Kd-Fab and BYbe**




Figure 38

**% Inhibition of Anti-IgM Induced CD40 on B-cells,
by Antigen 3 and Antigen 2 Specific Fab-Kd-Fab and BYbe**




Figure 39

**% Inhibition of Anti-IgM Induced CD86 on B-cells
by Antigen 3 and Antigen 2 Specific BYbe and Fab-Kd-Fab**




Figure 40




Figure 41

Percentage inhibition of CD71 expression on total B cells
stimulated with a goat anti human anti-IgM antibody




Figure 42




Figure 43

Percentage inhibition of CD27 expression on total B cells
stimulated with a goat anti human anti-IgM antibody




Figure 44




Figure 45

Percentage inhibition of CD86 expression on total B cells
stimulated with a goat anti human anti-IgM antibody

