



(12) **EUROPEAN PATENT SPECIFICATION**

(45) Date of publication and mention
of the grant of the patent:
18.10.2017 Bulletin 2017/42

(51) Int Cl.:
C12N 15/11 ^(2006.01) **C12Q 1/70** ^(2006.01)
C12N 15/867 ^(2006.01) **A61K 48/00** ^(2006.01)

(21) Application number: **14191299.8**

(22) Date of filing: **16.07.2011**

(54) **Retrovirus detection**

Retrovirusdetektion

Détection de rétrovirus

(84) Designated Contracting States:
**AL AT BE BG CH CY CZ DE DK EE ES FI FR GB
GR HR HU IE IS IT LI LT LU LV MC MK MT NL NO
PL PT RO RS SE SI SK SM TR**

(30) Priority: **16.07.2010 US 365297 P**
27.09.2010 US 386941 P
08.10.2010 US 391360 P

(43) Date of publication of application:
02.09.2015 Bulletin 2015/36

(62) Document number(s) of the earlier application(s) in
accordance with Art. 76 EPC:
11807616.5 / 2 580 325

(73) Proprietor: **Tocagen Inc.**
San Diego, CA 92109 (US)

(72) Inventors:
• **Jolly, Douglas, J.**
Encinitas, CA 92024 (US)
• **Perez, Omar**
San Diego, CA 92101 (US)
• **Lin, Amy**
San Diego, CA 92130 (US)

(74) Representative: **Daniels, Jeffrey Nicholas et al**
Page White & Farrer
Bedford House
John Street
London WC1N 2BF (GB)

(56) References cited:
WO-A2-2010/075414 WO-A2-2012/058637

- **URISMAN A ET AL:** "Identification of a novel Gammaretrovirus in prostate tumors of patients homozygous for R462Q RNASEL variant", PLOS PATHOGENS, PUBLIC LIBRARY OF SCIENCE, US; US NATIONAL LIBRARY OF MEDICINE (NLM), BETHESDA, MD, US, vol. 2, no. 3, 1 March 2006 (2006-03-01), pages E25-1, XP002669943, ISSN: 1553-7374, DOI: 10.1371/JOURNAL.PPAT.0020025
- **DATABASE Geneseq [Online] 22 May 2000** (2000-05-22), "Moloney murine leukaemia virus LTR PCR primer SEQ ID NO:16.", XP002739249, retrieved from EBI accession no. GSN:AAZ92907 Database accession no. AAZ92907 & US 6 015 694 A (DUBENSKY JR THOMAS W [US] ET AL) 18 January 2000 (2000-01-18)
- **DATABASE Geneseq [Online] 22 January 2009** (2009-01-22), "pGEM T easy-3'dLTR3 DNA amplifying PCR primer 3LTR-3, SEQ ID:9.", XP002739250, retrieved from EBI accession no. GSN:ATS37006 Database accession no. ATS37006 & KR 2008 0032562 A (VIROMED LTD [KR]) 15 April 2008 (2008-04-15)
- **DATABASE Geneseq [Online] 10 July 2008** (2008-07-10), "Viral microarray oligonucleotide sequence SEQ ID NO 11.", XP002739251, retrieved from EBI accession no. GSN:AXX99908 Database accession no. AXX99908 & WO 2007/130519 A2 (US GOV HEALTH & HUMAN SERV [US]; BAPTISTA CASSIO S [US]; WU XIAOLIN [U] 15 November 2007 (2007-11-15)

Note: Within nine months of the publication of the mention of the grant of the European patent in the European Patent Bulletin, any person may give notice to the European Patent Office of opposition to that patent, in accordance with the Implementing Regulations. Notice of opposition shall not be deemed to have been filed until the opposition fee has been paid. (Art. 99(1) European Patent Convention).

- DATABASE Geneseq [Online] 9 September 2004 (2004-09-09), "Oligonucleotide of the invention 3'-E/S/P/P-LINKER.", XP002739252, retrieved from EBI accession no. GSN:ADP67227 Database accession no. ADP67227 & WO 2004/052297 A2 (CYTOGENIX INC [US]; YIN CHEN [US]) 24 June 2004 (2004-06-24)
- OTTO ERLWEIN ET AL: "Failure to Detect the Novel Retrovirus XMRV in Chronic Fatigue Syndrome", PLOS ONE, vol. 5, no. 1, 6 January 2010 (2010-01-06) , page e8519, XP055072557, DOI: 10.1371/journal.pone.0008519
- VINCENT C LOMBARDI ET AL: "Detection of an infectious retrovirus, XMRV, in blood cells of patients with chronic fatigue syndrome", SCIENCE, AMERICAN ASSOCIATION FOR THE ADVANCEMENT OF SCIENCE, US; US NATIONAL LIBRARY OF MEDICINE (NLM), BETHESDA, MD, US, vol. 326, no. 5952, 23 October 2009 (2009-10-23), pages 585-589, XP002669944, ISSN: 1095-9203, DOI: 10.1126/SCIENCE.1179052 [retrieved on 2009-10-08]
- ROBERT SCHLABERG ET AL: "XMRV is present in malignant prostatic epithelium and is associated with prostate cancer, especially high-grade tumors", PROCEEDINGS OF THE NATIONAL ACADEMY OF SCIENCES, NATIONAL ACADEMY OF SCIENCES, US, vol. 106, no. 38, 22 September 2009 (2009-09-22), pages 16351-16356, XP007915265, ISSN: 0027-8424
- FISCHER N ET AL: "Prevalence of human gammaretrovirus XMRV in sporadic prostate cancer", JOURNAL OF CLINICAL VIROLOGY, ELSEVIER, AMSTERDAM, NL, vol. 43, no. 3, 1 November 2008 (2008-11-01), pages 277-283, XP025612273, ISSN: 1386-6532, DOI: 10.1016/J.JCV.2008.04.016 [retrieved on 2008-09-27]
- HOHN OLIVER ET AL: "Lack of evidence for xenotropic murine leukemia virus-related virus(XMRV) in German prostate cancer patients", RETROVIROLOGY, BIOMED CENTRAL LTD., LONDON, GB, vol. 6, no. 1, 16 October 2009 (2009-10-16), page 92, XP021063734, ISSN: 1742-4690, DOI: 10.1186/1742-4690-6-92
- B. DONG ET AL: "From the Cover: An infectious retrovirus susceptible to an IFN antiviral pathway from human prostate tumors", PROCEEDINGS OF THE NATIONAL ACADEMY OF SCIENCES, vol. 104, no. 5, 30 January 2007 (2007-01-30), pages 1655-1660, XP55015348, ISSN: 0027-8424, DOI: 10.1073/pnas.0610291104
- ARNOLD R S ET AL: "XMRV Infection in Patients With Prostate Cancer: Novel Serologic Assay and Correlation With PCR and FISH", UROLOGY, BELLE MEAD, NJ, US, vol. 75, no. 4, 1 April 2010 (2010-04-01), pages 755-761, XP026990712, ISSN: 0090-4295 [retrieved on 2010-04-01]

Description

TECHNICAL FIELD

[0001] The invention relates to methods and compositions useful for the detection of retroviruses in a subject or sample.

BACKGROUND

[0002] Xenotropic murine leukemia virus-related virus (XMRV) is a recently discovered human gammaretrovirus that resembles a xenotropic MLV, but that is distinguishable from xenotropic MLV in the sequence in its envelope (Urisman *et al.*, PLOS pathogens 2(3):e25, 2006; and Dong *et al.* PNAS 104:1655, 2007). All isolates so far examined are highly homologous to each other (>98% sequence identity) and allow the distinction from xenotropic MLV. The reason for this sequence conservation is not currently understood. The original infectious clone is called XMRV VP62 (GenBank accession no. EF185282).

[0003] XMRV was originally described in association with prostate cancer and further connections have been suggested (R. Schlaberg *et al.* PNAS 2009, doi_10.1073_pnas.0906922106), with 6-23% of prostate cancer patients testing positive. In addition, V.C.Lombardi *et al.* (Science 2009 doi10.1126/ science.11179052) showed a possible association with chronic fatigue syndrome, with 67% of patients testing positive, compared to 3.7% of normals. Overall estimates of the prevalence in the general population from investigators in the USA range from 2-4%. However several recent studies in Europe have failed to detect XMRV in similar frequencies or similar associations (Fischer *et al.* Journal of Clinical Virology 43: 277-283 2008; Hohn *et al.* Retrovirology 6:92 2009; FJM van Kuppeveld *et al.*, BMJ 340:c1018, 2010). Fischer *et al.* found 1 of 105 prostate cancer patients and 1 of 70 control subjects to be XMRV positive in a German study. Hohn *et al.* screened 589 prostate cancer patients in Germany without detecting a single positive. Van Kuppeveld *et al.* also failed to detect any DNA or RNA positives in 32 chronic fatigue patients or in 42 matched controls in Holland. Recently another paper from Fischer *et al.* (Emerg Infect Dis. 2010) showed about 10% positivity in RNA derived from sputum of 162 immunosuppressed patients and 2-3% positivity in sputum RNA from 168 normal patients in a German study. The assays did not appear to be different and no explanation was offered for the discrepancies.

[0004] The inconsistency of results calls into question the reliability of the current testing methods, in particular in DNA amplification. Following the lead of Lombardi *et al.* all investigators so far have used nested PCR using XMRV sequence based primers, followed by running the sample on a gel and looking for a visible band. This method is known to be variable in sensitivity and depend on the quality of nucleic acid samples. Detection of XMRV RNA has also been described mainly using the method of Dong *et al.* PNAS 2007, based on that of Urisman *et al.* 2006. In this assay RNA is prepared from tissue and/or blood, reverse transcribed to cDNA and the cDNA examined by QPCR with XMRV specific primers. As noted this led to inconsistent results (Enserink *et al.*, Science 329:18-19, 2010). Claims of various sensitivities have been made for such tests, but it is not possible to verify any of these and the assays appear to be incompletely characterized.

[0005] A PCR based diagnostic screening assay for XMRV in human blood has been recently developed (www.vip-dx.com), using nested PCR and gel detection of the amplification product (Lombardi *et al.*), with an estimated sensitivity for the nested DNA PCR around 600 proviral copies/test. In addition the report of Lombardi *et al.* do not show complete concordance of *gag* and *env* detection, with positives in *gag* and negatives for *env* observed in some subjects. This was attributed to variability in the assay. In all of the assays developed so far great care has been taken to use primers that will differentiate MLV from XMRV, so that only XMRV is detected. Therefore there is a great need for a reliable and validated assay for XMRV DNA and RNA in accessible samples from volunteers or patients in order to determine the real frequencies of positivity and whether there is linkage to disease. In addition a reliable blood screening assay is not available. Recent data suggest that detection of XMRV in many cases is caused by artifacts (Paprotka T., Science, 333, 97 - 101, 2011) or contamination with mouse DNA (Robinson MJ. *et al.*, Retrovirology, 7:108 doi:10.1186/1742-4690-7-108, 2010).

[0006] Furthermore, gene therapy vectors based upon MLV are being used including replication competent MLV-based vectors. For example, a replicating retrovirus based on amphotropic MLV and carrying an extra cytosine deaminase gene as a therapeutic agent for cancer including primary brain cancer leading to glioblastoma multiforme (GBM) (Tai *et al.*, Mol. Ther., 12:842-851 2005; http://oba.od.nih.gov/oba/RAC/meetings/Jun2009/976_Aghi.pdf; WO2010036986) having been used. An exemplary vector is being developed by TocaGen Inc. (San Diego, CA) and is referred to as Toca 511 (clinical trials.gov trial# NCT01156584). Subsequent to Toca 511 administration, patients are dosed with 5-fluorocytosine that is converted in situ to 5-fluorouracil, a potent anticancer compound. As the virus is generally only able to replicate in the tumor, this results in a very specific anti-cancer effect. In order to determine whether there is replication outside the tumor, for safety and/or for correlation with efficacy assays for detection of proviral DNA in the blood and MLV RNA in the plasma are needed. FDA currently requires follow-up on patients undergoing such investigational therapies with an integrating viral vector for 15 years post-treatment (Guidance for Industry - Supplemental Guidance

on Testing for Replication Competent Retrovirus in Retroviral Vector Based Gene Therapy Products and During Follow-up of Patients in Clinical Trials Using Retroviral Vectors: FDA Center for Biologics Evaluation and Research November 2006; <http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Cellular-andGeneTherapy/ucm072961.htm>). Assays are needed to accomplish this and a generally accepted marker for risk of disease and disease progression in viral diseases in general and retroviral diseases in particular is the levels over time of virus in the blood or in blood cells (Gurunathan S, Habib RE, et al. Vaccine. 2009;27:1997-2015; Low A., Okeoma CM. et al. Virology 2009; 385: 455-463). On the other hand, replication of the virus in the tumor may leak into the periphery and blood stream and so assays that monitor the appearance and levels of viral sequence in the blood as DNA or RNA can be used to determine whether there is an effective treatment and whether there is a need to modify the treatment protocol, for example to readminister the viral vector or to use adjuvants (such as steroids) that will facilitate the viral replication in the tumor.

[0007] WO 2010/075414 discloses the identification of Xenotropic murine leukaemia virus (XMRV) nucleic acid by polymerase chain reaction (PCR) analysis (e.g., real time PCR (RT/PCR); nested RT/PCR using Tth DNA polymerase and Hot start polymerase) and the uses thereof. Methods for the detection, and in particular early detection, of XMRV in RNA isolated from samples (e.g., urine samples; expressed prostate secretion (EPS)) of prostate cancer patients and normal individuals is disclosed.

[0008] Urisman et al., in "Identification of a novel Gammaretrovirus in prostate tumours of patients homozygous for R462Q RNASEL variant", PLOS Pathogens, vol. 2 (3), 2006, page E25-1, suggest that xenotropic MuLV-related viruses can produce an authentic human infection, and implicate RNase L activity in the prevention or clearance of infection in vivo.

SUMMARY

[0009] In a first aspect the present invention provides a composition defined in claim 1. In a second aspect the present invention provides a method defined in claim 2. In a third aspect the present invention provides a method defined in claim 4. In a fourth aspect the present invention provides a kit defined in claim 6. In a fifth aspect the present invention provides a method defined in claim 9. Further features of the invention are defined in the dependent claims.

[0010] The details of one or more embodiments of the disclosure are set forth in the accompanying drawings and the description below. Other features, objects, and advantages will be apparent from the description and drawings, and from the claims.

BRIEF DESCRIPTION OF DRAWINGS

[0011]

Figure 1 shows XMRV gag standard in pUC57-XMRV gag plasmid DNA. The insert corresponds to nucleotides 628 to 764 of the XMRV VP62 clone sequence (NC_007815). The derived sequence is synthesized by BioBasic Inc and inserted into pUC57 backbone at SmaI site between BamHI and ApaI sites.

Figure 2 shows XMRV env standard sequence in pET28b-XMRV env plasmid DNA. The insert corresponds to nucleotides 6252 to 6391 of the XMRV VP62 clone sequence (NC_007815). The derived sequence is synthesized by BioBasic Inc and inserted into pET28b+ backbone at EcoRV site between BssHI and HpaI sites.

Figure 3 shows the MLV pol and pol2 standard sequences in pAZ3-emd plasmid DNA which encodes an ecotropic Moloney MLV gag-pol, amphotropic env and IRES-GFPemd cassette downstream of the env (Logg et al. J.Virol. 75:6989-6998, 2001).

Figure 4 shows a comparison of the sequences of Toca 511 (the proviral form) used to treat GBM patients and an XMRV provirus (VP 62, NCBI Reference Sequence NC_007815.1), noting the overall homology of the LTR, gag, pol and envelope regions. Also shown are the regions of the pol1 and pol2 amplicons. The MLV pol1 and pol2 standard sequences in pAZ3-emd plasmid DNA.

Figure 5 shows a BLAST nucleic acid sequence comparison of the sequence of XMRV (VP62, NCBI ref. NC_007815.1, Sbjct) and MoMLV (NCBI ref. NC_001501.1, Query), showing sequences of 20 or more nucleotides that are exactly homologous (underlined/highlighted).

Figure 6 shows XMRV gag standard in pUC57-XMRV gag plasmid DNA amplification curves from 1E0 to 1E7.

Figure 7 shows XMRV env standard in pET28b-XMRV env plasmid DNA amplification curves from 1E0 to 1E7.

Figure 8A-C shows a) 1-Stage qPCR Protocol: pUC57-XMRV gag plasmid standard targeted with XMRV gag primer/probe set. pUC57-XMRV gag plasmid DNA in TE was targeted with XMRV gag primers and a 1-stage qPCR protocol was performed. The mean Ct and standard deviation was calculated; b) shows 1-Stage qPCR Protocol: Controls targeted with XMRV gag primer/probe set. 22Rv1 gDNA positive control, naïve human blood gDNA negative control and NTC were targeted with the XMRV gag primer/probe set and a 1-stage qPCR protocol was performed. The mean Ct, standard deviation and copies/reaction were calculated. 'ND' means 'non-detected'; c) shows 1-Stage

qPCR Protocol: Spiked human blood gDNA targeted with XMRV gag primer/probe set. Neat human blood gDNA was spiked with 8 log concentrations of pUC57-XMRV gag plasmid DNA (1E0 to 1E8 copies/reaction). The samples were targeted with the XMRV gag primer/probe set and a 1-stage qPCR protocol was performed. The mean Ct, standard deviation, copies/reaction and % recovery of the input copies/reaction were determined (the % recovery was determined by using the following equation: detected copies/reaction divided by the input copies/reaction times 100).

Figure 9A-C shows a) 1-Stage qPCR Protocol: pET28b-XMRV env plasmid standard targeted with XMRV env primer/probe set. pET28b-XMRV env plasmid DNA in TE was targeted with XMRV env primers and a 1-stage qPCR protocol was performed. The mean Ct and standard deviation was calculated; b) shows 1-Stage qPCR Protocol: Controls targeted with XMRV env primer/probe set. 22Rv1 gDNA positive control, naïve human blood gDNA negative control and NTC were targeted with the XMRV env primer/probe set and a 1-stage qPCR protocol was performed. The mean Ct, standard deviation and copies/reaction were calculated. 'ND' means 'non-detected'; c) shows 1-Stage qPCR Protocol: Spiked human blood gDNA targeted with XMRV env primer/probe set. Neat human blood gDNA was spiked with 8 log concentrations of pET28b-XMRV env plasmid DNA (1E0 to 1E8 copies/reaction). The samples were targeted with the XMRV env primer/probe set and a 1-stage qPCR protocol was performed. The mean Ct, standard deviation, copies/reaction and % recovery of the input copies/reaction were determined (the % recovery was determined by using the following equation: detected copies/reaction divided by the input copies/reaction times 100).

Figure 10A-C shows a) 1-Stage qPCR Protocol: pAZ3-emd pol2 plasmid standard targeted with XMRV pol2 primer/probe set. pAZ3-emd pol2 plasmid DNA in TE was targeted with XMRV pol2 primers and a 1-stage qPCR protocol was performed. The mean Ct and standard deviation was calculated; b) shows 1-Stage qPCR Protocol: Controls targeted with XMRV pol2 primer/probe set. 22Rv1 gDNA positive control, and naïve human blood gDNA negative control were targeted with the XMRV pol2 primer/probe set and a 1-stage qPCR protocol was performed. The mean Ct, standard deviation and copies/reaction were calculated. 'ND' means 'non-detected'; c) shows 1-Stage qPCR Protocol: Spiked human blood gDNA targeted with XMRV pol2 primer/probe set. Neat human blood gDNA was spiked with 8 log concentrations of pAZ3-emd pol2 plasmid DNA (1E0 to 1E8 copies/reaction). The samples were targeted with the XMRV pol2 primer/probe set and a 1-stage qPCR protocol was performed. The mean Ct, standard deviation, copies/reaction and % recovery of the input copies/reaction were determined (the % recovery was determined by using the following equation: detected copies/reaction divided by the input copies/reaction times 100).

Figure 11A-B shows a) 0-Stage vs. 1-Stage qPCR Protocols: pUC57 XMRV gag Standards. A 0-stage and a 1-stage qPCR protocol were performed targeting the pUC57 XMRV gag plasmid using XMRV gag primers. 'pUC57 XMRV gag' means the number of pUC57 XMRV gag copies spiked into a single qPCR reaction; b) shows 0-Stage vs. 1-Stage qPCR Protocols: pUC57 XMRV gag spike-ins into CA Human Blood gDNA. A 0-stage and a 1-stage qPCR protocol were performed targeting pUC57 XMRV gag spike-ins into CA human blood gDNA and using XMRV gag primers. 'pUC57 XMRV gag/001 gDNA' means the number of pUC57 XMRV gag copies spiked into donor 001 gDNA in a single qPCR reaction; '001' means 'donor #001'; 'ND' means 'non-detected'.

Figure 12A-B shows a) 0-Stage vs. 1-Stage qPCR Protocols: pET28b XMRV env Standards. A 0-stage and a 1-stage qPCR protocol were performed targeting the pET28b XMRV env plasmid using XMRV env primers. 'pET28b XMRV env' means the number of pET28b XMRV env copies spiked into a single qPCR reaction; b) shows 0-Stage vs. 1-Stage qPCR Protocols: pET28b XMRV env spike-ins into 001 Human Blood gDNA. A 0-stage and a 1-stage qPCR protocol were performed targeting pET28b XMRV env spike-ins into 001 human blood gDNA and using XMRV env primers. 'pET28b XMRV env /001 gDNA' means the number of pET28b XMRV env copies spiked into donor 001 gDNA in a single qPCR reaction; '001' means 'donor #001'; 'ND' means 'non-detected'.

Figure 13A-B shows a) 0-Stage vs. 1-Stage qPCR Protocols: pAZ3-emd pol2 standards. A 0-stage and a 1-stage qPCR protocol were performed targeting the pAZ3-emd pol2 plasmid using XMRV pol2 primers. 'pAZ3-emd pol2' means the number pAZ3-emd pol2 copies spiked into a single qPCR reaction; 'ND' means 'non-detected'; b) Shows 0-Stage vs. 1-Stage qPCR Protocols: pAZ3-emd pol2 spike-ins into 001 Human Blood gDNA. A 0-stage and a 1-stage qPCR protocol were performed targeting pAZ3-emd pol2 spike-ins into 001 human blood gDNA and using XMRV pol2 primers. 'pAZ3-emd pol2 /001 gDNA' means the number of pAZ3-emd pol2 copies spiked into donor 001 gDNA in a single qPCR reaction; '001' means 'donor #001'; 'ND' means 'non-detected'.

Figure 14 shows detection of MLV using MLV and ENV2 primer sets from formalin fixed paraffin embedded tissue (FFPE) infected with MLV. Paz3-emd spike in was added to either 100ng fresh tumor sample that was frozen or added to 100ng of a FFPE DNA tumor sample. qPCR was performed with the MLV and ENV2 primer sets.

Figure 15 shows detection of XMRV in whole blood by RTPCR using XMRV specific primer sets XMRV gag, XMRV pol2, XMRV env.

Figure 16A-B shows the results of monitoring patients over time with assays described herein for provirus DNA (MLVLTR primers and probes) in whole blood DNA, for viral RNA (by env RT-PCR) in the plasma, and for antiviral antibody responses in the plasma. These subjects (recurrent Glioblastoma multiforme (GBM) patients) were treated

by intracranial injection of 2.6×10^3 TU/g brain of T5.0002 amphotropic MLV retrovirus encoding a modified yeast cytosine deaminase (WO2010036986, WO2010045002) followed by 5-fluorocytosine treatment courses at approximately 130mg/kg/day. (A) patient 101; (B) patient 102.

DETAILED DESCRIPTION

[0012] Also, the use of "or" means "and/or" unless stated otherwise. Similarly, "comprise," "comprises," "comprising" "include," "includes," and "including" are interchangeable and not intended to be limiting.

[0013] It is to be further understood that where descriptions of various embodiments use the term "comprising," those skilled in the art would understand that in some specific instances, an embodiment can be alternatively described using language "consisting essentially of" or "consisting of."

[0014] Unless defined otherwise, all technical and scientific terms used herein have the same meaning as commonly understood to one of ordinary skill in the art to which this disclosure belongs. Although any methods and reagents similar or equivalent to those described herein can be used in the practice of the disclosed methods and compositions, the exemplary methods and materials are now described.

[0015] As used herein and in the appended claims, the singular forms "a," "and," and "the" include plural referents unless the context clearly dictates otherwise. Thus, for example, reference to "an oligonucleotide" includes a plurality of such oligonucleotides and reference to "the polynucleotide" includes reference to one or more polynucleotides known to those skilled in the art, and so forth.

[0016] The detection of XMRV or MLV related retroviruses by nucleic acid amplification techniques in human or animal tissues, blood or plasma/serum is of use for determining prostate cancer and risk thereof, chronic fatigue system and risk thereof, contamination of blood supply and tissue donation material and in following the status of subjects undergoing therapy with an MLV derived therapeutic virus comprising a heterologous genetic sequence such as, for example, an engineered retroviral replicatog virus based on amphotropic MLV (e.g., Toca 511). For example, the methods and compositions of the disclosure can be used to monitor therapy with an retroviral vector comprising sequences with substantial identity to MLV, in determining if recombination takes place between the therapeutic vector and XMRV or other MLV related natural infections, and for determining if a subject carries XMRV or another MLV related naturally occurring virus. Such assays are also useful for screening the blood supply to exclude subjects that are positive for XMRV or other MLV related retroviruses. Such assays also can be used to determine levels of MLV related virus over time, and provide information when it would be useful to start administering antiretroviral therapies that are also active against MLV such as, for example, AZT (Sakuma et al., Virology, 2009; Powell et al., J.Virol., 73:8813-8816, 1999; G.B.Beck-Engeser, PNAS, 2009). Such assays when used with histopathology samples can be used to determine the presence or absence of XMRV or other MLV related retroviruses in a patients stored sample or to determine the epidemiology of the XMRV or MLV related virus. Such assays can also be used to monitor patients to whom therapeutic vectors based on replicating MLV vectors have been administered. These measurements can be used to track the safety of the therapy over time (e.g., to 15 years and beyond) as high persistent levels (greater than 30,000, 100,000 or 300,000 copies/microgram) of MLV in genomic DNA or greater than 30,000 100,000 or 300,000 RNA copies/ml plasma) or increasing levels of these over time, can be used as a signal to more closely monitor for diseases that could be secondary to a therapy using an gene therapy vector comprising MLV or MLV-related sequences, such as leukemia or to start antiretroviral therapy. However, these measurements can also be used to judge the extent of replication of the MLV or MLV-related vector in a target tissue (i.e., efficacy or susceptibility to successful treatment) because of the possibility of "spill" into the circulatory system. Other uses of these assays for clinical monitoring will be apparent to those skilled in the art.

[0017] Engineered retroviral vectors that can be monitored include those set forth below:

RCR Vector - pAC-yCD2

tagttattaa tagtaatcaa ttacgggggtc attagttcat agcccatata tggagttccg
 cgttacataa cttacggtaa atggcccgcc tggctgaccg cccaacgacc cccgccatt
 gacgtcaata atgacgtatg ttcccatagt aacgccaata gggactttcc attgacgtca
 atgggtggag tatttacggg aaactgccc cttggcagta catcaagtgt atcatatgcc
 5 aagtacgccc cctattgacg tcaatgacgg taaatggccc gcctggcatt atgccagta
 catgacctta tgggactttc ctacttggca gtacatctac gtattagtca tcgctattac
 catggtgatg cggttttggc agtacatcaa tgggcgtgga tagcggtttg actcacgggg
 atttccaagt ctccacccca ttgacgtcaa tgggagtttg ttttggcacc aaaatcaacg
 ggactttcca aaatgtcgta acaactcgc cccattgacg caaatgggcg gtaggcgtgt
 10 acggtgggag gtctatataa gcagagctgg tttagtgaac cggcgccagt cctccgattg
 actgagtcgc ccgggtaccc gtgtatccaa taaacctct tgcagttgca tccgacttgt
 ggtctcgctg ttctttggga gggctctctc tgagtgattg actaccgctc agcgggggtc
 tttcatttgg gggctcgctc gggatcgga gacccctgcc cagggaccac cgaccacca
 ccgggaggtg agctggccag caacttatct gtgtctgtcc gattgtctag tgtctatgac
 15 tgattttatg cgcctgcgtc ggtactagtt agctaactag ctctgtatct ggcggacccg
 tgggtggaact gacgagttcg gaacaccg cgcgaacct gggagacgtc ccagggactt
 cgggggccgt ttttgtggcc cgacctgagt ccaaaaatcc cgatcgtttt ggactctttg
 gtgcaccccc cttagaggag ggatatgtgg ttctggtagg agacgagaac ctaaaacagt
 tcccgcctcc gtctgaattt ttgctttcgg tttgggaccg aagccgcgcc gcgcgtcttg
 20 tctgctgcag catcgttctg tgttgtctct gtctgactgt gtttctgtat ttgtctgaga
 atatgggcca gactgttacc actcccttaa gtttgacctt aggtcactgg aaagatgtcg
 agcggatcgc tcacaaccag tcggtagatg tcaagaagag acgttggggt accttctgct
 ctgcagaatg gccaaccttt aacgtcggat ggccgcgaga cggcaccttt aaccgagacc
 tcatcaccca ggttaagatc aaggtctttt cacctggccc gcatggacac ccagaccagg
 25 tcccctacat cgtgacctgg gaagccttgg cttttgaccc cctccctgg gtcaagccct
 ttgtacaccc taagcctccg cctcctcttc ctccatccgc cccgtctctc ccccttgaac
 ctccctgctt gacccgcct cgatcctccc tttatccage cctcactcct tctctaggcg
 ccaaacctaa acctcaagtt ctttctgaca gtggggggcc gctcatcgac ctacttacag
 aagaccccc gccttatagg gacccaagac caccctctc cgacagggac ggaaatgggtg
 30 gagaagcgac cctgcggga gaggcaccgg accctcccc aatggcatct cgcctacgtg
 ggagacggga gccccctgtg gccgactcca ctacctgca ggcattcccc ctccgcgcag
 gaggaacgg acagcttcaa tactggcgt tctcctctt tgacctttac aactggaaaa
 ataataaccc ttctttttct gaagatccag gtaaaactgac agctctgac gagtctgttc
 tcatcaccca tcagcccacc tgggacgact gtcagcagct gttggggact ctgctgaccg
 35 gagaagaaaa acaacgggtg ctcttagagg ctagaaaggc ggtgcggggc gatgatgggc
 gccccactca actgcccact gaagtcatg ccgcttttcc cctcgagcgc ccagactggg
 attacaccac ccaggcaggt aggaaccacc tagtccacta tcgccagttg ctccctagcg
 gtctccaaaa cgcgggcaga agccccacca atttggccaa ggtaaaagga ataacacaag
 ggcccaatga gtctccctcg gccttcctag agagacttaa ggaagcctat cgcaggtaca
 ctcccttatga ccctgaggac ccagggcaag aaactaatgt gtctatgtct ttcatttggc
 40 agtctgcccc agacattggg agaaagttag agaggttaga agatttaaaa aacaagacgc
 ttggagattt ggttagagag gcagaaaaga tctttaataa acgagaaacc ccggaagaaa
 gagaggaacg tatcaggaga gaaacagagg aaaaagaaga acgccgtagg acagaggatg

45

50

55

EP 2 913 402 B1

	agcagaaaga	gaaagaaaga	gatcgttagga	gacatagaga	gatgagcaag	ctattggcca
	ctgtcgtttag	tggacagaaa	caggatagac	agggaggaga	acgaaggagg	tcccaactcg
	atcgcgacca	gtgtgcctac	tgcaaagaaa	aggggactg	ggctaaagat	tgtcccaaga
5	aaccacgagg	acctcgggga	ccaagacccc	agacctccct	cctgacccta	gatgactagg
	gaggtcaggg	tcaggagccc	ccccctgaac	ccaggataac	cctcaaagtc	gggggggcaac
	ccgtcacctt	cctggtagat	actggggccc	aacactccgt	gctgacccaa	aatcctggac
	ccctaagtga	taagtctgcc	tgggtccaag	gggctactgg	aggaaagcgg	tatcgctgga
	ccacggatcg	caaagtacat	ctagctaccg	gtaaggtcac	ccactctttc	ctccatgtac
	cagactgtcc	ctatcctctg	ttaggaagag	atltgtctgac	taaactaaaa	gccccaatcc
10	acttttagagg	atcaggagcc	caggttatgg	gaccaatggg	gcagcccctg	caagtgttga
	ccctaaatat	agaagatgag	catcggctac	atgagacctc	aaaagagcca	gatgtttctc
	tagggctccac	atggctgtct	gatttttctc	aggcctgggc	ggaaaccggg	ggcatgggac
	tggcagttcg	ccaagctcct	ctgatcatac	ctctgaaagc	aacctctacc	cccgtgtcca
	taaaacaata	ccccatgtca	caagaagcca	gactggggat	caagccccac	atacagagac
	tgttggacca	gggaataact	gtaccctgac	agtccccctg	gaacacgccc	ctgctaccgg
15	ttaagaaacc	agggactatg	gattataggc	ctgtccagga	tctgagagaa	gtcaacaagc
	gggtggaaga	catccacccc	accgtgcccc	acccttacaa	cctcttgagc	gggctccac
	cgteccacca	gtggtacact	gtgcttgatt	taaaggatgc	ctttttctgc	ctgagactcc
	accccaccag	tcagcctctc	ttcgcccttg	agtggagaga	tccagagatg	ggaatctcag
	gacaattgac	ctggaccaga	ctccacacag	gtttcaaaaa	cagtcaccac	ctgtttgatg
20	aggcactgca	cagagaccta	gcagacttcc	ggatccagca	cccagacttg	atcctgtctac
	agtacgtgga	tgacttactg	ctggccgcca	cttctgagct	agactgccaa	caagggtactc
	gggccctggt	acaaacccta	gggaacctcg	ggtatcgggc	ctcggccaaag	aaagcccaaa
	tttgccagaa	acaggtcaag	tatctggggg	atcttctaaa	agagggtcag	agatggctga
	ctgaggccag	aaaagagact	gtgatggggc	agcctactcc	gaagaccctc	cgacaactaa
	gggagttcct	agggacggca	ggcttctgtc	gcctctggat	ccctgggttt	gcagaaatgg
25	cagccccctt	gtaccctctc	accaaaacgg	ggactctggt	taattggggc	ccagaccaac
	aaaaggccta	tcaagaaatc	aagcaagctc	ttctaactgc	cccagccctg	gggttgccag
	atltgactaa	gccctttgaa	ctctttgtcg	acgagaagca	gggctacgcc	aaagggtgtcc
	taacgcaaaa	actgggacct	tggcgtcggc	cggtggccta	cctgtccaaa	aagctagacc
	cagtagcagc	tgggtggccc	ccttgccctac	ggatggtagc	agccattgcc	gtactgacaa
30	aggatgcagg	caagctaacc	atgggacagc	cactagtcac	tctggccccc	catgcagtag
	aggcactagt	caaacaaccc	cccgaccgct	ggctttccaa	cgcccggtatg	actcactatc
	aggccttgct	tttgacacg	gaccgggtcc	agttcggacc	ggtggtagcc	ctgaaccggg
	ctacgtctgt	cccactgcct	gaggaagggc	tgcaacacaa	ctgccttgat	atcctggccg
	aagcccacgg	aacccgacct	gacctaacgg	accagccgct	cccagacgcc	gaccacacct
35	ggtagcacgga	tggaaagcgt	ctcttacaag	agggacagcg	taaggcgggg	gctgcgggtga
	ccaccgagac	cgaggtaatc	tgggctaaag	ccctgccagc	cgggacatcc	gctcagcggg
	ctgaactgat	agcactcacc	caggccctaa	agatggcaga	aggtaagaag	ctaaatgttt
	atactgatag	ccgttatgct	tttgctactg	cccatatcca	tggagaaata	tacagaaggc
	gtgggttgct	cacatcagaa	ggcaaagaga	tcaaaaataa	agacgagatc	ttggccctac
	taaaagccct	ctttctgccc	aaaagactta	gcataatcca	ttgtccagga	catcaaaagg
40	gacacagcgc	cgaggctaga	ggcaaccgga	tggctgacca	agcggcccga	aaggcagcca
	tcacagagac	tccagacacc	tctaccctcc	tcatagaaaa	ttcatcacc	tacacctcag
	aacatttttca	ttacacagtg	actgatataa	aggacctaac	caagttgggg	gccattttatg
	ataaaaacaaa	gaagtattgg	gtctaccaag	gaaaacctgt	gatgcctgac	cagtttactt
	ttgaattatt	agactttctt	catcagctga	ctcacctcag	cttctcaaaa	atgaaggctc
45	tcctagagag	aagccacagt	ccctactaca	tgtcgaaccg	ggatcgaaca	ctcaaaaata
	tactgagac	ctgcaaagct	tgtgcacaag	tcaacgccag	caagtctgcc	gttaaacagg
	gaactagggg	ccgcgggcat	cggcccgcca	ctcattggga	gatcgatttc	accgagataa
	agcccggtat	gtatggctat	aaatatcttc	tagtttttat	agataccttt	tctggctgga
	tagaagcctt	ccaaccaag	aaagaaaccg	ccaaggctcg	aaccaagaag	ctactagagg
	agatcttccc	caggttcggc	atgcctcagg	tattgggaac	tgacaatggg	cctgccttcg
50	tctccaaggt	gagtcagaca	gtggccgatc	tgttggggat	tgattggaaa	ttacattgtg
	catacagacc	ccaaagctca	ggccaggtag	aaagaatgaa	tagaaccatc	aaggagactt
	taactaaatt	aacgcttgca	actggctcta	gagactgggt	gctcctactc	cccttagccc
	tgtaccgagc	ccgcaacacg	ccgggcccc	atggcctcac	cccatatgag	actttatatg
	gggcaccccc	gccccctgtg	aacttccctg	accctgacat	gacaagagtt	actaacagcc
55	cctctctcca	agctcactta	caggctctct	acttagtcca	gcacgaagtc	tggagacctc
	tggcggcagc	ctaccaagaa	caactggacc	gaccgggtgg	acctcaccct	taccgagtcg
	gcgacacagt	gtgggtccgc	cgacaccaga	ctaagaacct	agaacctcgc	tggaaaggac

	cttacacagt	cctgctgacc	acccccacog	ccctcaaagt	agacggcatc	gcagcttgga
	tacacgccgc	ccacgtgaag	gctgccgacc	ccgggggtgg	accatcctct	agactgacat
	ggcgcggttca	acgctctcaa	aacccccctca	agataagatt	aaccctgtgga	agcccttaat
5	agtcattggga	gtcctgttag	gagtagggat	ggcagagagc	ccccatcagg	tctttaatgt
	aacctggaga	gtcaccaacc	tgatgactgg	gcgtaccgcc	aatgccacct	ccctcctggg
	aactgtacaa	gatgccttcc	caaaattata	ttttgatcta	tgtgatctgg	tcggagagga
	gtggggaccct	tcagaccagg	aaccgtatgt	cgggtatggc	tgcaagtacc	ccgcagggag
	acagcggacc	cggacttttg	acttttacgt	gtgccctggg	cataccgtaa	agtcgggggtg
	tgggggacca	ggagagggct	actgtggtaa	atgggggtgt	gaaaccaccg	gacaggctta
10	ctggaagccc	acatcatcgt	gggacctaat	ctcccttaag	cgcgtaaca	ccccctggga
	cacgggatgc	tctaaagttg	cctgtggccc	ctgctacgac	ctctccaaag	tatccaattc
	cttccaaggg	gctactcgag	ggggcagatg	caaccctcta	gtcctagaat	tcactgatgc
	aggaaaaaag	gctaactggg	acgggcccac	atcgtgggga	ctgagactgt	accggacagg
	aacagatcct	attaccatgt	tctccctgac	ccggcaggtc	cttaatgtgg	gaccccgagt
15	ccccataggg	cccaaccagc	tattaccoga	ccaaagactc	ccttcctcac	caatagagat
	tgtaccggct	ccacagccac	ctagccccct	caataccagt	tacccccctt	ccactaccag
	tacaccctca	acctccccct	caagtccaag	tgtcccacag	ccacccccag	gaactggaga
	tagactacta	gctctagtca	aaggagccta	tcaggcgctt	aacctcacca	atcccgacaa
	gacccaagaa	tggttggtgt	gcttagtgct	gggacctcct	tattacgaag	gagtagcggt
	cgtgggcact	tataccaatc	attccaccgc	tcgggccaac	tgtacggcca	cttcccaaca
20	taagcttacc	ctatctgaag	tgacaggaca	gggcctatgc	atgggggcag	tacctaaaac
	tcaccaggcc	ttatgtaaca	ccacccaaag	cgccggctca	ggatcctact	accttgcgagc
	acccgccgga	acaatgtggg	cttgccagcac	tggtattgact	ccttgcttgt	ccaccacggt
	gctcaatcta	accacagatt	atttgtgtatt	agttgaactc	tggcccagag	taattttacca
	ctcccccgat	tatatgtatg	gtcagcttga	acagcgtacc	aaatataaaa	gagagccagt
	atcattgacc	ctggcccttc	tactaggagg	attaaccatg	ggagggattg	cagctggaat
25	agggacgggg	accactgcct	taattaaaac	ccagcagttt	gagcagcttc	atgccgctat
	ccagacagac	ctcaacgaag	tcgaaaagtc	aattaccaac	ctagaaaagt	cactgacctc
	gttgtctgaa	gtagtcttac	agaaccgcag	aggcctagat	ttgctattcc	taaaggaggg
	aggtctctgc	gcagccctaa	aagaagaatg	ttgtttttat	gcagaccaca	cggggctagt
	gagagacagc	atggccaaat	taagagaaag	gcttaatcag	agacaaaaac	tatttgagac
30	aggccaagga	tggttcgaag	ggctgtttaa	tagatcccc	tggtttacca	ccttaatctc
	caccatcatg	ggacctctaa	tagtactctt	actgatctta	ctctttggac	cttgcatctt
	caatcgattg	gtccaatttg	ttaaagacag	gatctcagtg	gtccaggctc	tggttttgac
	tcagcaatat	caccagctaa	aacccataga	gtacgagcca	tgaacgcgtt	actggccgaa
	gccgcttgga	ataaggccgg	tgtgcgtttg	tctatatgtt	attttccacc	atattgccgt
35	cttttggtcaa	tgtgagggcc	cggaacctg	gccctgtctt	cttgacgagc	atctcaggg
	gtctttcccc	tctcgccaaa	ggaatgcaag	gtctgttgaa	tgtcgtgaag	gaagcagttc
	ctctggaagc	ttcttgaaag	caaacaacgt	ctgtagcgac	cctttgcagg	cagcggaacc
	ccccacctgg	cgacaggtgc	ctctgcggcc	aaaagccacg	tgtataagat	acacctgcaa
	aggcggcaca	acccagtgcc	cacgttgtga	gttgatagtg	tgtggaaaga	gtcaaatggc
	tctcctcaag	cgtattcaac	aaggggctga	aggatgcca	gaaggtaccc	cattgtatgg
40	gatctgatct	ggggcctcgg	tgcacatgct	ttacatgtgt	ttagtcgagg	ttaaaaaac
	gtctaggccc	cccgaaccac	ggggacgtgg	ttttcctttg	aaaaacacga	ttataaatgg
	tgaccggcgg	catggcctcc	aagtgggatc	aaaaggcat	ggatatcgct	tacgaggagg
	ccctgctggg	ctacaaggag	ggcggcgtgc	ctatcggcgg	ctgtctgatc	aacaacaagg
	acggcagtg	gctgggcagg	ggccacaaca	tgaggttcca	gaagggtcc	gccacctgc
45	acggcgagat	ctccacctg	gagaactgtg	gcaggctgga	gggcaagggtg	tacaaggaca
	ccacctgtga	caccacctg	tccccctgtg	acatgtgtac	cggcgctatc	atcatgtacg
	gcatccctag	gtgtgtgatc	ggcgagaacg	tgaacttcaa	gtccaagggtc	gagaagtacc
	tgcaaaccag	gggccacgag	gtggtggttg	ttgacgatga	gaggtgtaag	aagctgatga
	agcagttcat	cgacgagagg	cctcaggact	ggttcgagga	tatcggcgag	taagcggccg
	cagataaaat	aaaagatttt	atttagtctc	cagaaaaagg	ggggaatgaa	agacccacc
50	tgtaggtttg	gcaagctagc	ttaagtaacg	ccatttttgca	aggcatggaa	aaatacataa
	ctgagaatag	agaagttcag	atcaagggtc	ggaacagatg	gaacagctga	atatgggcca
	aacaggatat	ctgtggttaag	cagttcctgc	cccggctcag	ggccaagaac	agatggaaca
	gctgaatatg	ggccaaacag	gatatctgtg	gtaagcagtt	cctgccccgg	ctcagggcca
	agaacagatg	gtccccagat	gcggtccagc	cctcagcagt	ttctagagaa	ccatcagatg
55	tttccagggt	gccccaaagga	cctgaaatga	ccctgtgcct	tatttgaaact	aaccaatcag
	ttcgcttctc	gcttctgttc	gcgcgcttct	gctccccgag	ctcaataaaa	gagcccacaa
	ccctcactc	ggggcgccag	tcctccgatt	gactgagtcg	cccgggtacc	cgtgtatcca

EP 2 913 402 B1

ataaaccctc ttgcagttgc atccgacttg tgggtctcgct gttccttggg aggggtctcct
 ctgagtgatt gactaccggt cagcgggggt ctttcattac atgtgagcaa aaggccagca
 aaaggccagg aaccgtaaaa aggcgcggtt gctggcggtt tcccataggc tccgcccccc
 tgacgagcat cacaaaaatc gacgctcaag tcagaggtgg cgaaaccgga caggactata
 5 aagataccag gcgtttcccc ctggaagctc cctcgtgcgc tctcctgttc cgaccctgcc
 gcttaccgga tacctgtccg cctttctccc ttcggaagc gtggcgcttt ctcaatgctc
 acgctgtagg tatctcagtt cgggtgtaggt cgttcgctcc aagctgggct gtgtgcacga
 accccccgtt cagcccgacc gctgcgcctt atccggtaac tatcgtcttg agtccaaccc
 ggtaagacac gacttatcgc cactggcagc agccactggt aacaggatta gcagagcgag
 10 gtatgtaggc ggtgctacag agttcttgaa gtgggtggcct aactacggct acactagaag
 gacagtatct ggtatctgcy ctctgctgaa gccagttacc ttcggaaaaa gaggttggtag
 ctcttgatcc ggcaaacaaa ccaccgctgg tagcgggtgg ttttttggtt gcaagcagca
 gattacgcgc agaaaaaaag gatctcaaga agatcctttg atcttttcta cgggggtctga
 cgctcagtggt aacgaaaact cacgttaagg gatttttggtc atgagattat caaaaaggat
 15 cttcacctag atccttttaa attaaaaatg aagtttttaa tcaatctaaa gtatatatga
 gtaaacttgg tctgacagtt accaatgctt aatcagtgag gcacctatct cagcgatctg
 tctatttctg tcatccatag ttgcctgact ccccgctcgtg tagataacta cgatacggga
 gggcttacca tctggcccca gtgctgcaat gataccgcga gaccacgct caccggctcc
 agatttatca gcaataaacc agccagccgg aaggggccgag cgcagaagtg gtccctgcaac
 20 tttatccgcc tccatccagt ctattaattg ttgccgggaa gctagagtaa gtagttcgcc
 agttaatagt ttgcgcaacg ttgttgccat tgctgcaggc atcgtgggtg cacgctcgctc
 gtttggtatg gcttcattca gctccggttc ccaacgatca aggcgaggtta catgatcccc
 catgttggtgc aaaaaagcgg ttagctcctt cggctcctccg atcgttgtca gaagtaagtt
 ggccgcagtg ttatcactca tggttatggc agcactgcat aattctctta ctgtcatgcc
 25 atccgtaaga tgcttttctg tgactggtga gtactcaacc aagtcattct gagaatagtg
 tatgcggcga ccgagttgct cttgcccggc gtcaacacgg gataataccg cgccacatag
 cagaacttta aaagtgtca tcattggaaa acgttcttcg gggcgaaaac tctcaaggat
 cttaccgctg ttgagatcca gttcgatgta acccactcgt gcacccaact gatcttcagc
 atcttttact ttcaccagcg tttctgggtg agcaaaaaca ggaaggcaaa atgccgcaaa
 aaagggaata agggcgacac ggaaatggtt aatactcata ctcttccttt ttcaatatta
 30 ttgaagcatt tatcagggtt attgtctcat gagcggatac atatttgaat gtatttagaa
 aaataaacia ataggggttc cgcgcacatt tccccgaaaa gtgccacctg acgtctaaga
 aaccattatt atcatgacat taacctataa aaataggcgt atcacgaggc cctttcgtct
 tcaagaattc at

35 RCR Vector - pAC-yCD

40

45

50

55

EP 2 913 402 B1

	tagttattaa	tagtaatcaa	ttacgggggc	attagttcat	agcccatata	tggagttccg
	cgttacataa	cttacggtaa	atggcccgcc	tggctgaccg	cccaacgacc	cccggccatt
	gacgtcaata	atgacgtatg	ttcccatagt	aacgccaata	gggactttcc	attgacgtca
5	atgggtggag	tattttacgg	aaactgccc	cttggcagta	catcaagtgt	atcatatgcc
	aagtacgccc	cctattgacg	tcaatgacgg	taaatggccc	gcctggcatt	atgccagta
	catgacctta	tgggactttc	ctacttggca	gtacatctac	gtattagtca	tcgctattac
	catgggtgatg	cggttttggc	agtacatcaa	tgggcgtgga	tagcggtttg	actcacgggg
	atttccaagt	ctccacccca	ttgacgtcaa	tgggagtttg	ttttggcacc	aaaatcaacg
10	ggactttcca	aaatgtcgt	acaactccgc	cccattgacg	caaattggcg	gtaggcgtgt
	acggtgggag	gtctatataa	gcagagctgg	tttagtgaac	cggcgccagt	cctccgattg
	actgagtcgc	ccgggtaccc	gtgtatccaa	taaaccctct	tgcagttgca	tcgcacttgt
	ggtctcgctg	ttccttggga	gggtctcctc	tgagtgattg	actaccgctc	agcggggggtc
	tttcatttgg	gggctcgtcc	gggatcggga	gaccctgccc	cagggaccac	cgacccacca
15	ccgggaggta	agctggccag	caacttatct	gtgtctgtcc	gattgtctag	tgtctatgac
	tgattttatg	cgctgcgtc	ggtactagtt	agctaactag	ctctgtatct	ggcggaccgg
	tgggtggaact	gacgagttcg	gaacaccggg	ccgcaaccct	gggagacgtc	ccagggactt
	cgggggccgt	ttttgtggcc	cgacctgagt	ccaaaaatcc	cgatcgtttt	ggactctttg
	gtgcaccccc	cttagaggag	ggatatgtgg	ttctggtagg	agacgagaac	ctaaaacagt
20	tcccgcctcc	gtctgaattt	ttgctttcgg	tttgggaccg	aagccgcgcc	gcgcgtcttg
	tctgctgcag	catcgttctg	tgttgtctct	gtctgactgt	gtttctgtat	ttgtctgaga
	atatgggcca	gactgttacc	actcccttaa	gtttgacctt	aggtcactgg	aaagatgtcg
	agcggatcgc	tcacaaccag	tcggtagatg	tcaagaagag	acgttgggtt	accttctgct
	ctgcagaatg	gccaaccttt	aacgtcggat	ggccgcgaga	cggcaccttt	aaccgagacc
25	tcacaccca	ggttaagatc	aaggtctttt	cacctggccc	gcatggacac	ccagaccagg
	tcccctacat	cgtgacctgg	gaagccttgg	cttttgaccc	ccctccctgg	gtcaagccct

30

35

40

45

50

55

EP 2 913 402 B1

	ttgtacaccc	taagcctccg	cctcctcttc	ctccatccgc	cccgtctctc	ccccttgaac
	ctcctcgttc	gaccccgctt	cgatcctccc	tttatccagc	cctcactcct	tctctagggc
	ccaaacctaa	acctcaagtt	ctttctgaca	gtggggggcc	gctcatcgac	ctacttacag
5	aagaccccc	gccttatagg	gacccaagac	cacccccctc	cgacagggac	ggaaatggtg
	gagaagcgac	ccctgcggga	gaggcaccgg	acccctcccc	aatggcatct	cgcctacgtg
	ggagacggga	gccccctgtg	gccgactcca	ctacctcgca	ggcattcccc	ctccgcgcag
	gaggaaacgg	acagcttcaa	tactggccgt	tctcctcttc	tgacctttac	aactggaaaa
	ataataaccc	ttctttttct	gaagatccag	gtaaaactgac	agctctgata	gagtctgttc
	tcatcaccca	tcagcccacc	tgggacgact	gtcagcagct	gltggggact	ctgctgaccg
10	gagaagaaaa	acaacgggtg	ctcttagagg	ctagaaaggc	ggtgcggggc	gatgatgggc
	gccccactca	actgccaat	gaagtcgatg	ccgcttttcc	cctcgagcgc	ccagactggg
	attacaccac	ccaggcaggt	aggaaccacc	tagtccacta	tcgccagttg	ctcctagcgg
	gtctccaaaa	cgcgggcaga	agccccacca	atltggccaa	ggtaaaaggga	ataacacaag
	ggcccaatga	gtctccctcg	gccttcctag	agagacttaa	ggaagcctat	cgcaggtaca
	ctccttatga	ccctgaggac	ccagggcaag	aaactaatgt	gtctatgtct	ttcatttggc
15	agtctgcccc	agacattggg	agaaagttag	agaggttaga	agatttaaaa	aacaagcgc
	ttggagattt	ggttagagag	gcagaaaaga	tctttaataa	acgagaaacc	ccggaagaaa
	gagaggaacg	tatcaggaga	gaaacagagg	aaaaagaaga	acgccgtagg	acagaggatg
	agcagaaaaga	gaaagaaaga	gatcgttaga	gacatagaga	gatgagcaag	ctattggcca
	ctgtcgttag	tggacagaaa	caggatagac	agggaggaga	acgaaggagg	tcccaactcg
20	atcgcgacca	gtgtgcctac	tgcaaagaaa	aggggactg	ggctaaagat	tgtcccaaga
	aaccacgagg	acctcgggga	ccaagacccc	agacctccct	cctgacccta	gatgactagg
	gaggtcaggg	tcaggagccc	ccccctgaac	ccaggataac	cctcaaagtc	ggggggcaac
	ccgtcacctt	cctggtagat	actggggccc	aacactccgt	gctgacccaa	aatcctggac
	ccctaagtga	taagtctgcc	tgggtccaag	gggctactgg	aggaaagcgg	tatcgctgga
25	ccacggatcg	caaagtacat	ctagctaccg	gtaaggctac	ccactctttc	ctccatgtac
	cagactgtcc	ctatcctctg	ttaggaagag	atlttgctgac	taaaactaaa	gccccaatcc
	actttgaggg	atcaggagcc	caggttatgg	gaccaatggg	gcagcccctg	caagtgttga
	ccctaaatat	agaagatgag	catcggttac	atgagacctc	aaaagagcca	gatgtttctc
	tagggctccac	atggctgtct	gattttcctc	aggcctgggc	ggaaaccggg	ggcatgggac
	tggcagttcg	ccaagctcct	ctgatcatac	ctctgaaagc	aacctctacc	cccgtgtcca
30	taaaacaata	ccccatgtca	caagaagcca	gactggggat	caagccccac	atacagagac
	tgtttggacca	gggaatactg	gtaccctgcc	agtccccctg	gaacacgccc	ctgctacccg
	ttaagaaacc	agggactaat	gattataggg	ctgtccagga	tctgagagaa	gtcaacaagc
	gggtggaaga	catccacccc	accgtgcccc	acccttacia	cctcttgagc	gggtccccc
	cgtccaccac	gtggtacact	gtgcttgatt	ttaaaggatgc	ctttttctgc	ctgagactcc
35	accccaccag	tcagcctctc	ctcgcccttg	agtgagagaga	tccagagatg	ggaatctcag
	gacaattgac	ctggaccaga	ctccacacagg	gtttcaaaaa	cagtcccacc	ctgtttgatg
	aggcactgca	cagagacctc	gcagacttcc	ggatccagca	cccagacttg	atcctgctac
	agtacgtgga	tgacttactg	ctggccgcca	cttctgagct	agactgccaa	caagggtactc
	gggcccctgtt	acaaacccta	gggaacctcg	ggtatcgggc	ctcggccaa	aaagcccaaa
	tttgccagaa	acaggtcaag	tatctggggg	atcttctaaa	agagggtcag	agatggctga
40	ctgaggccag	aaaagagact	gtgatggggc	agcctactcc	gaagaccocct	cgacaactaa
	gggagttcct	agggacggca	ggcttctgtc	gcctctggat	ccctgggttt	gcagaaatgg
	cagccccctt	gtaccctctc	accaaaacgg	ggactctgtt	taattggggc	ccagaccaac
	aaaaggccta	tcaagaaatc	aagcaagctc	ttctaactgc	cccagccctg	gggttgccag
	atlttgactaa	gccctttgaa	ctctttgtcg	acgagaagca	gggtacgccc	aaaggtgtcc
45	taacgcaaaa	actgggacct	tggcgctcgg	cgggtggccta	cctgtccaaa	aagctagacc
	cagtagcagc	tgggtggccc	ccttgccctac	ggatggtagc	agccattgcc	gtactgacaa
	aggatgcagg	caagctaacc	atgggacagc	cactagtcat	tctggccccc	catgcagtag
	aggcactagt	caaacaaccc	cccgaccgct	ggctttccaa	cgcccgatg	actcactatc
	aggccttgct	tttgacacg	gaccgggtcc	agttcggacc	ggtggtagcc	ctgaaccggy
	ctacgctgct	cccactgcct	gaggaagggc	tgcaacacaa	ctgccttgat	atcctggccg
50	aagcccacgg	aacccgaccc	gacctaacgg	accagccgct	cccagacgcc	gaccacacct
	ggtacacgga	tggaaagcag	ctcttacaag	agggacagcg	taaggcgggga	gctgcgggtga
	ccaccgagac	cgaggtaatc	tgggctaagc	ccctgccagc	cgggacatcc	gctcagcggg
	ctgaactgat	agcactcacc	caggccctaa	agatggcaga	aggtaaagaa	ctaaatgttt
	atactgatag	ccgttatgct	tttgctactg	cccatatcca	tggagaaata	tacagaaggc
55	gtgggttgct	cacatcagaa	ggcaaagaga	tcaaaaaata	agacgagatc	ttggccctac
	taaaagccct	ctttctgccc	aaaagactta	gcataatcca	ttgtccagga	catcaaaagg
	gacacagcgc	cgaggctaga	ggcaaccgga	tggctgacca	agcggcccga	aaggcagcca

EP 2 913 402 B1

	tacacagagac	tccagacacc	tctaccctcc	tcatagaaaa	ttcatcacc	tacacctcag
	aacatttttca	ttacacagt	actgatataa	aggacctaac	caagttgggg	gccattttatg
	ataaaaacaaa	gaagtatttg	gtctaccaag	gaaaacctgt	gatgcctgac	cagttttactt
5	ttgaattatt	agactttctt	catcagctga	ctcacctcag	cttctcaaaa	atgaaggctc
	tcttagagag	aagccacagt	ccctactaca	tgctgaaccg	ggatcgaaca	ctcaaaaata
	tcactgagac	ctgcaaagct	tgtgcacaag	tcaacgccag	caagtctgcc	gttaaacagg
	gaactaggg	ccgcgggcat	cggcccggca	ctcattggga	gatcgatttc	accgagataa
	agcccggatt	gtatggctat	aaatatcttc	tagtttttat	agataccttt	tctggctgga
	tagaagcctt	cccaaccaag	aaagaaaccg	ccaaggctcg	aaccaagaag	ctactagagg
10	agatcttccc	caggttcggc	atgcctcagg	tattgggaac	tgacaatggg	cctgccttcg
	tctccaaggt	gagtcagaca	gtggccgatc	tggtggggat	tgattggaaa	ttacattgtg
	catacagacc	ccaaagctca	ggccaggtag	aaagaatgaa	tagaaccatc	aaggagactt
	taactaaatt	aacgcttgca	actggctcta	gagactgggt	gctcctactc	cccttagccc
	tgtaccgagc	ccgcaacacg	cggggccccc	atggcctcac	cccatatgag	atcttatatg
	gggcaccccc	gccccttgta	aacttccctg	accctgacat	gacaagagtt	actaacagcc
15	cctctctcca	agctcactta	caggctctct	acttagtcca	gcacgaagtc	tgagagcctc
	tgggcggcagc	ctaccaagaa	caactggacc	gaccgggtgt	acctcaccct	taccgagtcg
	gcgacacagt	gtgggtccgc	cgacaccaga	ctaagaacct	agaacctcgc	tggaaaggac
	cttacacagt	cctgctgacc	acccccaccg	ccctcaaagt	agacggcatc	gcagcttgga
	tacacgccgc	ccacgtgaag	gctgccgacc	ccgggggtgg	accatcctct	agactgacat
20	ggcgcggttca	acgctctcaa	aacccccctc	agataagatt	aaccctggga	agcccttaat
	agtcattggga	gtcctgttag	gagtagggat	ggcagagagc	ccccatcagg	tctttaatgt
	aacctggaga	gtcaccaacc	tgatgactgg	gcgtaccgcc	aatgccacct	ccctcctggg
	aactgtacaa	gatgccttcc	caaaattata	ttttgatcta	tgtgatctgg	tggagagga
	gtgggaccct	tcagaccagg	aaccgtatgt	cgggtatggc	tgcaagtacc	ccgcaggggag
	acagcggacc	cggacttttg	acttttacgt	gtgccctggg	cataccgtaa	agtcgggggtg
25	tgggggacca	ggagaggctt	actgtggtaa	atgggggtgt	gaaaccaccg	gacaggctta
	ctggaagccc	acatcatcgt	gggacctaat	ctcccttaag	cgcggttaaca	ccccctggga
	cacgggatgc	tctaaagttg	cctgtggccc	ctgctacgac	ctctccaaag	tatccaattc
	cttccaagg	gctactcgag	ggggcagatg	caaccctcta	gtcctagaat	tactgatgc
	agggaaaaaag	gctaactggg	acggggccca	atcgtgggga	ctgagactgt	accggacagg
30	aacagatcct	attaccatgt	tctccctgac	ccggcaggtc	cttaatgtgg	gaccccgagt
	ccccatagg	cccaaccag	tattacccga	ccaaagactc	ccttcctcac	caatagagat
	tgtaccggct	ccacagccac	ctagccccct	caataccagt	tacccccctt	ccactaccag
	tacacctcta	acctccccta	caagtccaag	tgtcccacag	ccacccccag	gaactggaga
	tagactacta	gctctagtca	aaggagccta	tcaggcgctt	aacctcacca	atcccagaca
35	gacccaagaa	tggttggtgt	gcttagtgct	gggacctcct	tattacgaag	gagtagcggt
	cgtgggcact	tataccaatc	attccaccgc	tccggccaac	tgtacggcca	cttcccaaca
	taagcttacc	ctatctgaag	tgacaggaca	gggcctatgc	atggggggcag	tacctaaaac
	tcaccaggcc	ttatgtaaca	ccacccaaag	cgccggctca	ggatcctact	accttgcagc
	acccgccgga	acaatgtggg	cttgcagcac	tggtattgact	ccctgcttgt	ccaccacggg
40	gctcaatcta	accacagatt	attgtgtatt	agttgaactc	tggccagag	taatttacca
	ctccccgat	tatatgtatg	gtcagcttga	acagcgtacc	aaatataaaa	gagagccagt
	atcattgacc	ctggcccttc	tactaggagg	attaaccatg	ggagggattg	cagctggaat
	agggacgggg	accactgcct	taattaaaac	ccagcagttt	gagcagcttc	atgccgctat
	ccagacagac	ctcaacgaag	tcgaaaagtc	aattaccaac	ctagaaaagt	cactgacctc
	gttgtctgaa	gtagtcctac	agaaccgcag	aggcctagat	ttgctattcc	taaaggagg
45	aggtctctgc	gcagccctaa	aagaagaatg	ttgtttttat	gcagaccaca	cggggctagt
	gagagacagc	atggccaaat	taagagaaag	gcttaatcag	agacaaaaac	tatttgagac
	aggccaagga	tggttcgaag	ggctgtttta	tagatcccc	tggtttacca	ccttaatctc
	caccatcatg	ggacctctaa	tagtactctt	actgatctta	ctctttggac	cttgcattct
	caatcgattg	gtccaatttg	ttaaagacag	gatctcagtg	gtccaggctc	tggttttgac
50	tcagcaatat	caccagctaa	aacccataga	gtacgagcca	tgaacgcgtt	actggccgaa
	gocgcttgga	ataaggccgg	tgtgcgtttg	tctatatgtt	attttccacc	atattgccgt
	cttttggcaa	tgtgagggcc	cggaaacctg	gccctgtctt	cttgacgagc	attcctaggg
	gtctttcccc	tctcgccaaa	ggaatgcaag	gtctgttgaa	tgtcgtgaag	gaagcagttc
	ctctggaagc	ttcttgaaga	caaacaacgt	ctgtagcgac	cctttgcagg	cagcggaacc
	ccccacctgg	cgacaggtgc	ctctgcggcc	aaaagccacg	tgtataagat	acacctgcaa
55	aggcggcaca	acccagtgct	cacgttgtga	gttgatagtg	tgtggaaaga	gtcaaatggc
	tctcctcaag	cgtattcaac	aaggggctga	aggatgcccc	gaaggtaccc	cattgtatgg
	gatctgatct	ggggcctcgg	tgacatgct	ttacatgtgt	ttagtcgagg	ttaaaaaac

EP 2 913 402 B1

	gtctaggccc	cccgaaccac	ggggacgtgg	ttttcctttg	aaaaacacga	ttataaatgg
	tgacaggggg	aatggcaagc	aagtgggatc	agaaggggat	ggacattgcc	tatgaggagg
	cggccttagg	ttacaaagag	ggtggtgttc	ctattggcgg	atgtcttata	aataacaaag
5	acggaagtgt	tctcggtcgt	ggtcacaaca	tgagatttca	aaagggatcc	gccacactac
	atggtgagat	ctccactttg	gaaaactgtg	ggagattaga	gggcaaagtg	tacaaagata
	ccactttgta	tacgacgctg	tctccatgcg	acatgtgtac	aggtgccatc	atcatgtatg
	gtattccacg	ctgtgttgct	ggtgagaacg	ttaatttcaa	aagtaagggc	gagaaatatt
	tacaaactag	aggtcacgag	gttggtgttg	ttgacgatga	gaggtgtaaa	aagatcatga
	aacaatttat	cgatgaaaga	cctcaggatt	ggtttgaaga	tattggtgag	tagggcgccg
10	cagataaaat	aaaagatttt	atttagtctc	cagaaaaagg	ggggaatgaa	agaccccacc
	tgtaggtttg	gcaagctagc	ttaagtaacg	ccattttgca	aggcatggaa	aaatacataa
	ctgagaatag	agaagttcag	atcaaggtca	ggaacagatg	gaacagctga	atatggggcca
	aacaggatat	ctgtggtaag	cagttcctgc	ccgggtcag	ggccaagaac	agatggaaca
	gctgaatatg	ggccaaacag	gatatactgtg	gtaagcagtt	cctgcccccg	ctcaggggcca
15	agaacagatg	gtccccagat	gcgggtccagc	cctcagcagt	ttctagagaa	ccatcagatg
	tttccagggt	gccccaaagg	cctgaaatga	ccctgtgcct	tatttgaact	aaccaatcag
	ttcgcttctc	gcttctgttc	gcgcgcttct	gctccccgag	ctcaataaaa	gagcccacaa
	cccctcactc	ggggcgccag	tcctccgatt	gactgagtcg	cccgggtacc	cgtgtatcca
	ataaaccctc	ttgcagttgc	atccgacttg	tggtctcgtc	gttccttggg	aggggtctcct
20	ctgagtgatt	gactaccggt	cagcgggggt	ctttcattac	atgtgagcaa	aaggccagca
	aaaggccagg	aaccgtaaaa	aggccgcgtt	gctggcggtt	ttccataggc	tcggcccccc
	tgacgagcat	cacaaaaatc	gacgctcaag	tcagaggtgg	cgaaccccg	caggactata
	aagataccag	gcgtttcccc	ctggaagctc	cctcgtgcgc	tctcctgttc	cgaccctgcc
	gcttaccgga	tacctgtccg	cctttctccc	ttcgggaagc	gtggcgcttt	ctcaatgctc
25	acgctgtagg	tatctcagtt	cgggtgtagg	cgttcgctcc	aagctgggct	gtgtgcacga
	acccccggtt	cagcccgacc	gctgcgcctt	atccggtaac	tatcgtcttg	agtccaaccc
	ggtaagacac	gacttatcgc	cactggcagc	agccactggt	aacaggatta	gcagagcgag
	gtatgtaggc	ggtgctacag	agttcctgaa	gtggtggcct	aactacggct	acactagaag
	gacagtattt	ggtatctgcg	ctctgctgaa	gccagttacc	ttcggaaaaa	gagttggtag
30	ctcttgatcc	ggcaaacaaa	ccaccgctgg	tagcgggtgg	ttttttgttt	gcaagcagca
	gattacgcgc	agaaaaaaag	gatctcaaga	agatcctttg	atcttttcta	cggggtctga
	cgctcagtg	aacgaaaact	cacgttaagg	gatttttggtc	atgagattat	caaaaaggat
	cttcacctag	atccttttaa	attaaaaatg	aagtttttaa	tcaatctaaa	gtatatatga
	gtaaacttgg	tctgacagtt	accaatgctt	aatcagtgag	gcacctatct	cagcgatctg
	tctatttcgt	tcatccatag	ttgcctgact	ccccgtcgtg	tagataacta	cgatacggga
35	gggcttacca	tctggcccca	gtgctgcaat	gataccgcga	gaccacgct	caccggctcc
	agatttatca	gcaataaacc	agccagccgg	aagggccgag	cgcagaagtg	gtcctgcaac
	tttatccgcc	tccatccagt	ctattaattg	ttgccgggaa	gctagagtaa	gtagttcgcc
	agttaatagt	ttgcgcaacg	ttgttgccat	tgctgcaggc	atcgtgggtg	cacgctcgct
	gtttggtag	gcttcattca	gctccgggtc	ccaacgatca	aggcgagtta	catgatcccc
40	catgttgtgc	aaaaaagcgg	ttagctcctt	cggctcctcg	atcgttgtca	gaagtaagtt
	ggcgcagtg	ttatcactca	tggttatggc	agcactgcat	aattctctta	ctgtcatgcc
	atccgtaaga	tgcttttctg	tgactgggtga	gtactcaacc	aagtcattct	gagaatagtg
	tatgcggcga	ccgagttgct	cctgcccggc	gtcaacacgg	gataataacc	cgccacatag
	cagaacttta	aaagtgtca	tcattggaaa	acgttcttcg	gggcgaaaac	tctcaaggat
45	cttaccgctg	ttgagatcca	gttcgatgta	accactcgt	gcacccaact	gatcttcagc
	atcttttact	ttcaccagcg	tttctgggtg	agcaaaaaa	ggaaggcaaa	atgccgcaaa
	aaagggaata	agggcgacac	ggaaatgttg	aataactata	ctcttccttt	ttcaatatta
	ttgaagcatt	tatcagggtt	attgtctcat	gagcggatag	atatttgaat	gtatttagaa
	aaataaacia	ataggggttc	cgcgcacatt	tccccgaaaa	gtgccacctg	acgtctaaga
50	aaccattatt	atcatgacat	taacctataa	aaataggcgt	atcacgaggc	cctttcgtct
	tcaagaattc	at				

RCR Vector - pACE-CD

EP 2 913 402 B1

tagttattaa tagtaatcaa ttacgggggtc attagttcat agcccatata tggagttccg
cgttacataa cttacggtaa atggcccgcc tggctgaccg cccaacgacc cccgcccatt
gacgtcaata atgacgtatg ttcccatagt aacgccaata gggactttcc attgacgtca
atgggtggag tatttacggt aaactgcca cttggcagta catcaagtgt atcatatgcc
aagtacgcc cctattgacg tcaatgacgg taaatggccc gcctggcatt atgcccagta
catgacctta tgggactttc ctacttggca gtacatctac gtattagtca tcgctattac
catggtgatg cggtttttggc agtacatcaa tgggcgtgga tagcggtttg actcacgggg

5

10

15

20

25

30

35

40

45

50

55

EP 2 913 402 B1

	atttccaagt	ctccacccca	ttgacgtcaa	tgggagtttg	ttttggcacc	aaaatcaacg
	ggactttcca	aaatgtcgtg	acaactccgc	cccattgacg	caaagtggcg	gtaggcgtgt
	acgggtgggag	gtctatataa	gcagagctgg	tttagtgaac	cggcgccagt	cctccgattg
5	actgagtcgc	ccgggtaccc	gtgtatccaa	taaaccctct	tgcagttgca	tccgacttgt
	ggctctcgctg	ttccttggga	gggtctcctc	tgagtgattg	actaccgctc	agcggggggtc
	tttcattttg	gggctcgtcc	gggatcgga	gaccctgccc	cagggaccac	cgaccacca
	ccgggaggtg	agctggccag	caacttatct	gtgtctgtcc	gattgtctag	tgtctatgac
	tgattttatg	cgctgcgtc	ggtactagtt	agctaactag	ctctgtatct	ggcggaccgc
	tgggtggaact	gacgagttcg	gaacaccccg	ccgcaaccct	gggagacgtc	ccagggactt
10	cgggggccgt	ttttgtggcc	cgacctgagt	ccaaaaatcc	cgatcgtttt	ggactctttg
	gtgcaccccc	cttagaggag	ggatatgtgg	ttctggtagg	agacgagaac	ctaaaacagt
	tcccgctcc	gtctgaattt	ttgctttcgg	tttgggaccg	aagccgcgcc	gcgcgtcttg
	tctgctgcag	catcgttctg	tgttgtctct	gtctgactgt	gtttctgtat	ttgtctgaga
	atatgggcca	gactgttacc	actcccttaa	gtttgacctt	aggtcactgg	aaagatgtcg
	agcggatcgc	tcacaaccag	tcggtagatg	tcaagaagag	acgttggggt	accttctgct
15	ctgcagaatg	gccaaccttt	aacgtcggat	ggccgcgaga	cggcaccttt	aaccagacc
	tcatacccca	ggttaagatc	aaggctcttt	cacctggccc	gcatggacac	ccagagagg
	tcccctacat	cgtgacctgg	gaagccttgg	cttttgacct	ccctccctgg	gtcaagccct
	ttgtacaccc	taagcctccg	cctcctcttc	ctccatccgc	cccgtctctc	ccccttgaac
	ctcctcgttc	gaccccgct	cgatcctccc	tttatccagc	cctcactcct	tctctaggcg
20	ccaaacctaa	acctcaagtt	ctttctgaca	gtggggggcc	gctcatcgac	ctacttacag
	aagaccccc	gccttatagg	gacccaagac	caccccttcc	cgacagggac	ggaaatggtg
	gagaagcgac	ccctgcggga	gaggcaccgg	acccctcccc	aatggcatct	cgctacgtg
	ggagacggga	gccccctgtg	gccgactcca	ctacctcgca	ggcattcccc	ctccgcgcag
	gaggaaacgg	acagcttcaa	tactggcgtg	tctcctcttc	tgacctttac	aactggaaaa
25	ataataaccc	ttctttttct	gaagatccag	gtaaaactgac	agctctgac	gagctgttct
	tcatacccca	tcagcccacc	tgggacgact	gtcagcagct	gttggggact	ctgctgaccg
	gagaagaaaa	acaacgggtg	ctcttagagg	ctagaaaggc	ggtgcggggc	gatgatgggc
	gccccactca	actgccaat	gaagtcgatg	ccgcttttcc	cctcgagcgc	ccagactggg
	attacaccac	ccaggcaggt	aggaaccacc	tagtccacta	tcgccagttg	ctcctagcgg
	gtctccaaaa	cgcgggcaga	agccccacca	atltggccaa	ggtaaaagga	ataacacaag
30	ggcccaatga	gtctccctcg	gccttcctag	agagacttaa	ggaagcctat	cgcaggtaca
	ctccttatga	ccctgaggac	ccagggcaag	aaactaatgt	gtctatgtct	ttcatttggc
	agtctgcccc	agacattggg	agaaaagttag	agaggttaga	agatttaaaa	aacaagacgc
	ttggagatlt	ggttagagag	gcagaaaaga	tctttaataa	acgagaaacc	ccggaagaaa
	gagaggaacg	tatcaggaga	gaaacagagg	aaaaagaaga	acgccgtagg	acagaggtg
35	agcagaaaga	gaaagaaaaga	gatcgtagga	gacatagaga	gatgagcaag	gatgtttcca
	ctgtcgttag	tggacagaaa	caggatagac	agggaggaga	acgaaggagg	tcccaactcg
	atcgcgacca	gtgtgcctac	tgcaaagaaa	aggggcaactg	ggctaagat	tgtcccaaga
	aaccacgagg	acctcgggga	ccaagacccc	agacctccct	cctgacccta	gatgactagg
	gaggtcaggg	tcaggagccc	ccccctgaac	ccaggataac	cctcaaagtc	ggggggcaac
	ccgtcacctt	cctggttagat	actggggccc	aacactccgt	gctgacccaa	aatcctggac
40	ccctaagtga	taagtctgcc	tgggtccaag	gggctactgg	aggaaagcgg	tatcgctgga
	ccacggatcg	caaagtacat	ctagctaccg	gtaaggctac	ccactctttc	ctccatgtac
	cagactgtcc	ctatcctctg	ttaggaagag	atltgtctgac	taaactaaaa	gccccaatcc
	actttgaggg	atcaggagcc	caggttatgg	gaccaatggg	gcagcccctg	caagtgttga
	ccctaaatat	agaagatgag	catcggtac	atgagacctc	aaaagagcca	gatgtttctc
45	tagggctccac	atggctgtct	gattttcctc	aggcctgggc	ggaaaccggg	ggcatgggac
	tggcagttcg	ccaagctcct	ctgatcatac	ctctgaaagc	aacctctacc	cccgtgtcca
	taaaacaata	ccccatgtca	caagaagcca	gactggggat	caagccccac	atacagagac
	tgtttggacca	gggaataactg	gtaccctgcc	agtccccctg	gaacacgccc	ctgctacccg
	ttaagaaacc	agggactaat	gattataggc	ctgtccagga	tctgagagaa	gtcaacaagc
	gggtggaaga	catccacccc	accgtgccca	acccttacia	cctcttgagc	gggctcccac
50	cgteccacca	gtggtacact	gtgcttgatt	taaaggatgc	ctttttctgc	ctgagactcc
	accccaccag	tcagcctctc	ttcgcccttg	agtggagaga	tccagagatg	ggaatctcag
	gacaattgac	ctggaccaga	ctcccacagg	gtttcaaaaa	cagtcaccac	ctggttgatg
	aggcacctga	agggactaat	gcagacttcc	ggatccagca	cccagacttg	atcctgtac
	agtacgtgga	tgacttactg	ctggccgcca	cttctgagct	agactgccaa	caaggctactc
55	gggcccgtgt	acaaacccta	gggaacctcg	ggatatgggc	ctcggccaa	aaagcccaaa
	tttgccagaa	acaggtcaag	tatctggggg	atcttctaaa	agagggtcag	agatggctga
	ctgaggccag	aaaagagact	gtgatggggc	agcctactcc	gaagaccctc	cgacaactaa

EP 2 913 402 B1

	gggagttcct	agggacggca	ggcttctgtc	gcctctggat	ccctgggttt	gcagaaatgg
	cagccccctt	gtaccctctc	acaaaaacgg	ggactctgtt	taattggggc	ccagaccaac
	aaaaggccta	tcaagaaatc	aagcaagctc	ttctaactgc	cccagccctg	gggttgccag
5	atttgactaa	gccctttgaa	ctctttgtcg	acgagaagca	gggtacgcc	aaaggtgtcc
	taacgcaaaa	actgggacct	tggcgctcggc	cgggtggccta	cctgtccaaa	aagctagacc
	cagtagcagc	tgggtggccc	ccttgccctac	ggatggtagc	agccattgcc	gtactgacaa
	aggatgcagg	caagctaacc	atgggacagc	cactagtcat	tctggccccc	catgcagtag
	aggcactagt	caaacaaccc	cccgaccgct	ggctttccaa	cgcccgatg	actcactatc
	aggccttgct	tttgacacg	gaccgggtcc	agttcggacc	ggtggtagcc	ctgaacccgg
10	ctacgtctgt	cccactgcct	gaggaagggc	tgcaacacaa	ctgccttgat	atcctggccg
	aagcccacgg	aacccgaccc	gacctaacgg	accagccgct	cccagacgcc	gaccacacct
	ggtacacgga	tggaaagcagt	ctcttacaag	agggacagcg	taaggcggga	gctgcgggtga
	ccaccgagac	cgaggtaatc	tgggctaaag	ccctgccagc	cgggacatcc	gctcagcggg
	ctgaactgat	agcactcacc	caggccctaa	agatggcaga	aggtaagaag	ctaaatgttt
	atactgatag	ccgttatgct	tttgctactg	cccatatcca	tggagaaata	tacagaaggc
15	gtgggttgct	cacatcagaa	ggcaaagaga	tcaaaaataa	agacgagatc	ttggccctac
	taaaagccct	ctttctgccc	aaaagactta	gcataatcca	ttgtccaggga	catcaaaagg
	gacacagcgc	cgaggctaga	ggcaaccgga	tggctgacca	agcggcccga	aaggcagcca
	tcacagagac	tccagacacc	tctaccctcc	tcatagaaaa	ttcatcacc	tacacctcag
	aacattttca	ttacacagtg	actgatataa	aggacctaac	caagttgggg	gccatttatg
20	ataaaacaaa	gaagtattgg	gtctaccaag	gaaaacctgt	gatgcctgac	cagtttactt
	ttgaattatt	agactttctt	catcagctga	ctcacctcag	cttctcaaaa	atgaaggctc
	tcctagagag	aagccacagt	ccctactaca	tgctgaaccg	ggatcgaaca	ctcaaaaata
	tactgagac	ctgcaaagct	tgtgcacaag	tcaacgccag	caagtctgcc	gttaaacagg
	gaactagggt	ccgcgggcat	cggcccggca	ctcattggga	gatcgatttc	accgagataa
	agcccggatt	gtatggctat	aaatatcttc	tagtttttat	agataccttt	tctggctgga
25	tagaagcctt	cccaaccaag	aaagaaaccg	ccaaggtcgt	aaccaagaag	ctactagagg
	agatcttccc	cagggttcggc	atgcctcagg	tattgggaac	tgacaatggg	cctgccttcg
	tctccaaggt	gagtcagaca	gtggccgcatc	tgttggggat	tgattggaaa	ttacattgtg
	catacagacc	ccaaagctca	ggccaggtag	aaagaatgaa	tagaaccatc	aaggagactt
	taactaaatt	aacgcttgca	actggctcta	gagactgggt	gctcctactc	cccttagccc
30	tgtaccgagc	ccgcaacacg	ccgggcccc	atggcctcac	cccatatgag	atccttatatg
	gggcaccccc	gccccttgta	aacttccctg	accctgacat	gacaagagtt	actaacagcc
	cctctctcca	agctcactta	caggtctctc	acttagtcca	gcacgaagtc	tggagacctc
	tggcggcagc	ctaccaagaa	caactggacc	gaccggtggt	acctcaccct	taccgagtcg
	gcgacacagt	gtgggtccgc	cgacaccaga	ctaagaacct	agaacctcgc	tggaaaggac
35	cttacacagt	cctgctgacc	acccccacgc	ccctcaaaag	agacgcatc	cagcttgga
	tacacgccgc	ccacgtgaag	gctgccgacc	ccgggggtgg	accatcctct	agactgacat
	ggcgcggttca	acgctctcaa	aacccccctca	agataagatt	aacccgtgga	agcccttaat
	agtcattggga	gtcctgttag	gagtagggat	ggcagagagc	ccccatcagg	tctttaatgt
	aacctggaga	gtcaccaacc	tgatgactgg	gcgtaccgcc	aatgccacct	ccctcctggg
	aactgtacaa	gatgccttcc	caaaattata	ttttgatcta	tgtgatctgg	tggagagga
40	gtgggacctt	tcagaccagg	aaccgtatgt	cgggtatggc	tgcaagtacc	ccgcagggag
	acagcggacc	cggacttttg	acttttacgt	gtgccctggg	cataccgtaa	agtcgggggtg
	tgggggacca	ggagagggt	actgtggtaa	atgggggtgt	gaaaccaccg	gacaggctta
	ctggaagccc	acatcatcgt	gggacctaat	ctcccttaag	cgcgtaaca	ccccctggga
	cacgggatgc	tctaaaagttg	cctgtggccc	ctgctacgac	ctctccaaa	tatccaattc
45	cttccaaggg	gctactcgag	ggggcagatg	caaccctcta	gtcctagaat	tcactgatgc
	aggaaaaaag	gctaactggg	acgggcccac	atcgtgggga	ctgagactgt	accggacagg
	aacagatcct	attaccatgt	tctccctgac	ccggcaggtc	cttaatgtgg	gaccccgagt
	ccccataggg	cccaaccag	tattaccoga	ccaaagactc	ccttcctcac	caatagagat
	tgtaccggct	ccacagccac	ctagccccct	caataccagt	tacccccctt	ccactaccag
	tacaccctca	acctccccca	caagtccaag	tgtcccacag	ccacccccag	gaactggaga
50	tagactacta	gctctagtca	aaggagccta	tcaggcgctt	aacctcacca	atcccgacaa
	gacccaagaa	tgttggtgtg	gcttagtgct	gggacctcct	tattacgaag	gagtgcgggt
	cgtgggcact	tataccaatc	attccaccgc	tccggccaac	tgtacggcca	cttcccaaca
	taagctttacc	ctatctgaag	tgacaggaca	gggcctatgc	atgggggcca	taacctaaac
	tcaccaggcc	ttatgtaaca	ccacccaaag	cgccggctca	ggatcctact	accttgagc
55	acccgcggga	acaatgtggg	cttgagcac	tggattgact	ccctgcttgt	ccaccacggg
	gctcaatcta	accacagatt	attgtgtatt	agttgaactc	tggcccagag	taatttacca
	ctcccccgat	tatatgtatg	gtcagcttga	acagcgtacc	aaatataaaa	gagagccagt

EP 2 913 402 B1

	atcattgacc	ctggcccttc	tactaggagg	attaaccatg	ggagggattg	cagctggaat
	agggacgggg	accactgcct	taattaaaac	ccagcagttt	gagcagcttc	atgccgctat
	ccagacagac	ctcaacgaag	tcgaaaagtc	aattaccaac	ctagaaaagt	cactgacctc
	gttgtctgaa	gtagtcctac	agaaccgcag	aggcctagat	ttgctattcc	taaaggaggg
5	aggctctctgc	gcagccctaa	aagaagaatg	ttgtttttat	gcagaccaca	cggggctagt
	gagagacagc	atggccaaat	taagagaaag	gcttaatcag	agacaaaaac	tatttgagac
	aggccaagga	tggttcgaag	ggctgtttaa	tagatcccc	tggtttacca	ccttaatctc
	caccatcatg	ggacctctaa	tagtactctt	actgatctta	ctctttggac	cttgcaattct
	caatcgatta	gtccaatttg	ttaaagacag	gatatcagtg	gtccaggctc	tagttttgac
10	tcaacaatat	caccagctga	agcctataga	gtacgagcca	tgacgtacgt	tactggccga
	agccgcttgg	aataaggccg	gtgtgcgttt	gtctatatgt	tattttccac	catattgccg
	tcttttggca	atgtgagggc	ccggaaacct	ggccctgtct	tcttgacgag	cattcctagg
	ggtctttccc	ctctcgccaa	aggaatgcaa	ggtctgttga	atgtcgtgaa	ggaagcagtt
	cctctggaag	cttcttgaag	acaaacaacg	tctgtagcga	ccctttgcag	gcagcggaac
	ccccacctg	gcgacaggtg	cctctgcggc	caaaagccac	gtgtataaga	tacacctgca
15	aaggcggcac	aaccccagtg	ccacgttgtg	agttggatag	ttgtgaaag	agtc aaatgg
	ctctcctcaa	gcgtattcaa	caaggggctg	aaggatgcc	agaaggtagc	ccattgtatg
	ggatctgata	tggggcctcg	gtgcacatgc	tttacctgtg	tttagtcgag	gttaaaaaaa
	cgtctaggcc	ccccgaacca	cggggacgtg	gttttctctt	gaaaaacacg	ataataccat
	ggtgacaggg	ggaatggcaa	gcaagtggga	tcagaagggt	atggacattg	cctatgagga
20	ggcggcctta	ggttacaaag	aggggtgggt	tcctattggc	ggatgtctta	tcaataacaa
	agacggaagt	gttctcgttc	gtggtcacaa	catgagattt	caaaagggat	ccgccacact
	acatggtgag	atctccactt	tggaaaactg	tgggagatta	gagggcaaa	tgtacaaaga
	taccactttg	tatacgacgc	tgtctccatg	cgacatgtgt	acaggtgcca	tcacatgta
	tggtattcca	cgctgtgttg	tcggtgagaa	cgtaaatttc	aaaagtaagg	gcgagaaata
	tttacaact	agaggtcacg	aggttgttgt	tgttgacgat	gagaggtgta	aaaagatcat
25	gaaacaattt	atcgatgaaa	gacctcagga	ttggtttgaa	gatattgggtg	agtaggcggc
	cgcgccatag	ataaaaataaa	agattttatt	tagtctccag	aaaaaggggg	gaatgaaaga
	ccccacctgt	aggtttggca	agctagctta	agtaacgcca	ttttgcaagg	catggaaaaa
	tacataactg	agaatagaga	agttcagatc	aaggtcagga	acagatggaa	cagctgaata
	tgggccaac	aggatatctg	tggtaagcag	ttcctgcccc	ggctcagggc	caagaacaga
30	tggaaacagct	gaatatgggc	caaacaggat	atctgtggta	agcagttcct	gccccggctc
	agggccaaaga	acagatggtc	cccagatgcg	gtccagccct	cagcagtttc	tagagaacca
	tcagatgttt	ccagggtgcc	ccaaggacct	gaaatgacct	tgtgccttgt	ttaaactaac
	caatcagttt	gcttctcgct	tctgttcgcg	cgcttctgct	ccccgagctc	aataaaaagag
	cccacaacct	ctcactcggg	gcgccagtc	tcggattgac	tgagtcgccc	gggtaccgct
35	gtatccaata	aacctcttgg	cagttgcata	cgacttgtgg	tctcgtgtgt	tcttgggagg
	gtctcctctg	agtgattgac	taccgcgtcag	cgggggtctt	tcatttgggg	gctcgtccgg
	gatcgggaga	ccctgcccc	gggaccaccg	accaccacc	gggaggtaa	ctggctgcct
	cgcgcggttt	ggtgatgacg	gtgaaaacct	ctgacatgtg	agcaaaaggc	cagcaaaagg
	ccaggaaccg	taaaaaggcc	gcgttgctgg	cgtttttcca	taggtccgc	ccccctgacg
	agcatcacaa	aaatcgacgc	tcaagtcaga	ggtggcgaaa	cccgacagga	ctataaagat
40	accaggcggt	tccccctgga	agctccctcg	tgcgctctcc	tgttccgacc	ctgccgctta
	ccggataacct	gtccgccttt	ctcccttcgg	gaagcgtggc	gctttctcaa	tgtctacgct
	gtaggatatct	cagttcggtg	taggtcggtc	gtcctaagct	gggctgtgtg	cacgaacccc
	cggttcagcc	cgaccgctgc	gccttatccg	gtaactatcg	tcttgagtc	aacctggtaa
	gacacgactt	atcgccactg	gcagcagcca	ctggtaacag	gattagcaga	gcgaggtatg
45	taggcgggtg	tacagagttc	ttgaagtgg	ggcctaacta	cggctacact	agaaggacag
	tatttgggat	ctgcgctctg	ctgaagccag	ttaccttcgg	aaaaagagtt	ggtagctctt
	gatccggcaa	acaaaccacc	gctggtagcg	gtgggttttt	tgtttgcaag	cagcagatta
	cgcgcgagaaa	aaaaggatct	caagaagatc	ctttgatctt	ttctacgggg	tctgacgctc
	agtggaacga	aaactcacgt	taagggtatt	tggtcatgag	attatcaaaa	aggatcttca
	cctagatcct	tttaaattaa	aaatgaagtt	ttaaatcaat	ctaaagtata	tatgagtaaa
50	cttgggtctga	cagttacca	tgtttaatca	gtgaggcacc	tatctcagcg	atctgtctat
	ttcgttcctc	catagttgcc	tgactccccg	tcgtgtagat	aactacgata	cgggagggct
	taccatctgg	ccccagtgct	gcaatgatac	cgcgagaccc	acgctcaccg	gctccagatt
	tatcagcaat	aaaccagcca	gcccgaagg	ccgagcgacg	aagtggctct	gcaactttat
	ccgcctccat	ccagttctatt	aattgttgcc	gggaagctag	agtaagtagt	tcgccagtta
55	atagtttgcg	caacgttggt	gccattgctg	caggcatcgt	ggtgtcacgc	tcgtcgtttg
	gtatggcttc	attcagctcc	ggttcccaac	gatcaaggcg	agttacatga	tccccatgt
	tgtgcaaaaa	agcggttagc	tccttcggtc	ctccgatcgt	tgtcagaagt	aagttggccg

EP 2 913 402 B1

5 cagtgttatac actcatgggtt atggcagcac tgcataattc tcttactgtc atgccatccg
taagatgctt ttctgtgact ggtgagtact caaccaagtc attctgagaa tagtgtatgc
ggcgaccgag ttgctcttgc ccggcgtcaa cacgggataa taccgcgcca catagcagaa
ctttaaaagt gctcatcatt ggaaaacgtt cttcggggcg aaaactctca aggatcttac
cgctggtgag atccagttcg atgtaaccca ctcggtgcacc caactgatct tcagcatctt
ttactttcac cagcgtttct gggtgagcaa aaacaggaag gcaaaatgcc gcaaaaaagg
gaataagggc gacacggaaa tgttgaatac tcatactctt cttttttcaa tattattgaa
gcatttatca gggttattgt ctcatgagcg gatacatatt tgaatgtatt tagaaaaata
10 aacaaatagg ggttccgcgc acatttcccc gaaaagtgcc acctgacgtc taagaaacca
ttattatcat gacattaacc tataaaaaata ggcgtatcac gaggcccttt cgtcttcaag
aattcat

RCR Vector - pAC-yCD

15

20

25

30

35

40

45

50

55

tagttattaa tagtaatcaa ttacgggggc attagttcat agcccatata tggagttccg
 cgttacataa cttacggtaa atggcccgcc tggctgaccg cccaacgacc cccgccatt
 gacgtcaata atgacgtatg ttcccatagt aacgccataa gggactttcc attgacgtca
 atgggtggag tattttacggt aaactgccc a ttggcagta catcaagtgt atcatatgcc
 5 aagtacgccc cctattgacg tcaatgacgg taaatggccc gcctggcatt atgccagta
 catgacctta tgggactttc ctacttggca gtacatctac gtattagtca tcgctattac
 catgggtgatg cggttttggc agtacatcaa tgggcgtgga tagcggtttg actcacgggg
 atttccaagt ctccacccca ttgacgtcaa tgggagtttg ttttggcacc aaaatcaacg
 ggactttcca aaatgtcgta acaactccgc ccattgacg caaatgggcg gtaggcgtgt
 10 acggtgggag gtctatataa gcagagctgg tttagtgaac cggcgccagt cctccgattg
 actgagtcgc ccgggtaccc gtgtatccaa taaacctct tgcagttgca tccgacttgt
 ggtctcgctg ttccctggga gggctctctc tgagtgattg actaccgctc agcgggggctc
 tttcatttgg gggctcgctc gggatcggga gacctctgcc cagggaccac cgaccacca
 ccgggaggt aagctggccag caacttatct gtgtctgtcc gattgtctag tgtctatgac
 15 tgattttatg cgctgcgctc ggtactagtt agctaactag ctctgtatct ggcggaccgg
 tgggtggaatg gacgagttcg gaacaccgg ccgcaaccct gggagacgct ccagggactt
 cggggggcgt ttttgtggcc cgacctgagt ccaaaaatcc cgatcgtttt ggactccttg
 gtgcaccccc cttagaggag ggatatgtgg ttctggtagg agacgagaac ctaaaacagt
 tccgcctcc gtctgaattt ttgtcttcgg tttgggaccg aagccgcgcc gcgcgtcttg
 20 tctgctgcag catcgttctg tgttgtctct gtctgactgt gtttctgtat ttgtctgaaa
 atatgggcca gactgttacc actcccttaa gtttgacctt aggtcactgg aaagatgtcg
 agcggatcgc tcacaaccag tcggtagatg tcaagaagag acgttgggtt accttctgct
 ctgcagaatg gccaaccttt aacgtcggat ggccgcgaga cggcaccttt aaccgagacc
 tcatcaccca ggttaagatc aaggtctttt cacctggccc gcatggacac ccagaccagg
 25 tcccctacat cgtgacctgg gaagccttgg cttttgacct cctccctgg gtcaagccct
 ttgtacaccc taagcctccg cctcctcttc ctccatccgc cccgtctctc ccccttgaac
 ctccctgctt gaccccgctt cgatcctccc tttatccagc cctcactcct tctctaggcg
 ccaaacctaa acctcaagtt ctttctgaca gtgggggggc gctcatcgac ctacttacag
 aagacccccc gccttatagg gacccaagac caccctcttc cgacagggac ggaaatgggtg
 30 gagaagcgac ccctgcggga gaggcaccgg accctcccc aatggcatct cgcctacgtg
 ggagacggga gcccctgtg gccgactcca ctacctcgca ggcattcccc ctccgcgag
 gaggaaacgg acagcttcaa tactggcgt tctcctcttc tgacctttac aactggaaaa
 ataataaccc ttctttttct gaagatccag gtaaaactgac agctctgac gatctgtctc
 tcatcaccca tcagcccacc tgggacgact gtcagcagct gttggggact ctgctgaccg
 35 gagaagaaaa acaacgggtg ctcttagagg ctagaaaggc ggtgcggggc gatgatgggc
 gccccactca actgcccact gaagtcatg ccgcttttcc cctcgagcgc ccagactggg
 attacaccac ccaggcaggt aggaaccacc tagtccacta tcgccagttg ctccctagcg
 gtctccaaaa cgcgggcaga agccccacca atttggccaa ggtaaaagga ataacacaag
 ggccaatga gtctccctcg gccttcttag agagacttaa ggaagcctat cgcaggtaca
 40 ctcttatga ccctgaggac ccagggaag aaactaatgt gtctatgtct ttcatttggc
 agtctgcccc agacattggg agaaagttag agaggtaga agatttaaaa aacaagacgc
 ttggagattt ggttagagag gcagaaaaga tctttaataa acgagaaacc ccggaagaaa
 gagaggaacg tatcaggaga gaaacagagg aaaaagaaga acgccgtagg acagaggatg
 agcagaaaaga gaaagaaaaga gatcgtagga gacatagaga gatgagcaag ctattggcca
 45 ctgtcgttag tggacagaaa caggatagac agggaggaga acgaaggagg tcccaactcg
 atcgcgacca gtgtgcctac tgcaaagaaa aggggcaact ggctaagat tgtcccaaga
 aaccacgagg acctcgggga ccaagacccc agacctcct cctgacccta gatgactagg
 gaggtcaggg tcaggagccc cccctgaac ccaggataac cctcaaagtc ggggggcaac

50

55

EP 2 913 402 B1

	cogtcacctt	cctggtagat	actggggccc	aacactccgt	gctgacccaa	aatcctggac
	ccctaagtga	taagtctgcc	tgggtccaag	gggctactgg	aggaaagcgg	tatcgctgga
	ccacggatcg	caaagtacat	ctagctaccg	gtaaggtcac	ccactctttc	ctccatgtac
	cagactgtcc	ctatcctctg	ttaggaagag	atltgtctgac	taaactaaaa	gccccaatcc
5	acttttaggg	atcaggagcc	caggttatgg	gaccaatggg	gcagcccctg	caagtgttga
	ccctaataat	agaagatgag	tatcggctac	atgagacctc	aaaagagcca	gatgtttctc
	taggggtccac	atggctgtct	gatttttctc	aggcctgggc	ggaaaccggg	ggcatgggac
	tggcagttcg	ccaagctcct	ctgatcatac	ctctgaaagc	aacctctacc	cccgtgtcca
	taaaacaata	ccccatgtca	caagaagcca	gactggggat	caagccccac	atacagagac
10	tgtttggacca	gggaatactg	gtaccctgcc	agtccccctg	gaacacgccc	ctgctacccg
	ttaagaaacc	agggactaat	gattataggc	ctgtccagga	tctgagagaa	gtcaacaagc
	gggtggaaga	catccacccc	accgtgccca	acccttacia	cctcttgagc	gggctcccac
	cgteccacca	gtggtacact	gtgcttgatt	taaaggatgc	ctttttctgc	ctgagactcc
	accccaccag	tcagcctctc	ttcgcccttg	agtggagaga	tccagagatg	ggaatctcag
	gacaattgac	ctggaccaga	ctcccacag	gtttcaaaaa	cagtcccacc	ctggttgatg
15	aggcactgca	cagagactta	gcagacttcc	ggatccagca	cccagacttg	atcctgttac
	agtagctgga	tgacttactg	ctggccgcca	cttctgagct	agactgccaa	caaggctactc
	gggcccgtgt	acaaacccta	gggaacctcg	ggatcgggc	ctcggccaag	aaagcccaaa
	tttgccagaa	acaggtcaag	tatctggggg	atcttctaaa	agagggtcag	agatggctga
	ctgaggccag	aaaagagact	gtgatggggc	agcctactcc	gaagaccctc	cgacaactaa
20	gggagttcct	agggacggca	ggcttctgtc	gcctctggat	ccctgggttt	gcagaaatgg
	cagccccctt	gtaccctctc	acaaaaacgg	ggactctgtt	taattggggc	ccagaccaac
	aaaaggccta	tcaagaaatc	aagcaagctc	ttctaactgc	cccagccctg	gggttgccag
	atltgactaa	gccctttgaa	ctctttgtcg	acgagaagca	gggctacgcc	aaagggtgtcc
	taacgcaaaa	actgggacct	tggcgctcgg	cgggtggccta	cctgtccaaa	aagctagacc
	cagtagcagc	tgggtggccc	ccttgccctac	ggatggtagc	agccattgcc	gtactgacaa
25	aggatgcagg	caagctaacc	atgggacagc	cactagtcac	tctggccccc	catgcagtag
	aggcactagt	caaacaaccc	cccgaccgct	ggctttccaa	cgcccgatg	actcactatc
	aggccttgct	tttgacacg	gaccgggtcc	agttcggacc	ggtggtagcc	ctgaaccggg
	ctacgtctgt	cccactgcct	gaggaagggc	tgcaacacaa	ctgccttgat	atcctggccg
	aagcccacgg	aacccgaccc	gacctaacgg	accagccgct	cccagacgcc	gaccacacct
30	ggtacacgga	tggaagcagt	ctcttacaag	agggacagcg	taaggcggga	gctgcgggtga
	ccaccgagac	cgaggtaatc	tgggctaaag	ccctgccagc	cgggacatcc	gctcagcggg
	ctgaactgat	agcactcacc	caggccctaa	agatggcaga	aggtaagaag	ctaaatgttt
	atactgatag	ccgttatgct	tttgctactg	cccatatcca	tggagaaata	tacagaaggc
	gtgggttgct	cacatcagaa	ggcaaagaga	tcaaaaataa	agacgagatc	ttggccctac
35	taaaagccct	ctttctgccc	aaaagactta	gcataatcca	ttgtccagga	ctacaaaagg
	gacacagcgc	cgaggctaga	ggcaaccgga	tggctgacca	agcggcccga	aaggcagcca
	tcacagagac	tccagacacc	tctaccctcc	tcatagaaaa	ttcatcacc	tacacctcag
	aacattttca	ttacacagtg	actgatataa	aggacctaac	caagttgggg	gccatttatg
	ataaaaacaa	gaagtattgg	gtctaccaag	gaaaacctgt	gatgcctgac	cagtttactt
	ttgaattatt	agactttctt	catcagctga	ctcacctcag	cttctcaaaa	atgaaggctc
40	tcctagagag	aagccacagt	ccctactaca	tgctgaaccg	ggatcgaaca	ctcaaaaata
	tactgagac	ctgcaaagct	tgtgcacaag	tcaacgccag	caagtctgcc	gttaaacagg
	gaactagggt	ccgcgggcat	cggcccgcca	ctcattggga	gatcgatttc	accgagataa
	agcccggatt	gtatggctat	aaatatcttc	tagtttttat	agataccttt	tctggctgga
	tagaagccct	cccaaccaag	aaagaaaccg	ccaaggtcgt	aaccaagaag	ctactagagg
45	agatcttccc	caggttcggc	atgcctcagg	tattgggaac	tgacaatggg	cctgccttcg
	tctccaagg	gagtcagaca	gtggccgata	tgttggggat	tgattggaaa	ttacattgtg
	catacagacc	ccaaagctca	ggccaggtag	aaagaatgaa	tagaaccatc	aaggagactt
	taactaaatt	aacgcttgca	actggctcta	gagactgggt	gctcctactc	cccttagccc
	tgtaccgagc	ccgcaacacg	ccgggcccc	atggcctcac	cccatatgag	atcttatatg
	gggcaccccc	gccccttgta	aacttccctg	accctgacat	gacaagagtt	actaacagcc
50	cctctctcca	agctcactta	caggtctctc	acttagtcca	gcacgaagtc	tggagacctc
	tggcggcagc	ctaccaagaa	caactggacc	gaccgggtgt	acctcaccct	taccgagtgc
	gcgacacagt	gtgggtccgc	cgacaccaga	ctaagaacct	agaacctcgc	tggaaaggac
	cttacacagt	cctgctgacc	acccccacgc	ccctcaaagt	agacggtgac	gagacttgga
	tacacgccgc	ccacgtgaag	gctgccgacc	ccgggggtgg	accatcctct	acactgacat
55	ggcgcggttca	acgctctcaa	aacccccctc	agataagatt	aacccgtgga	agcccttaat
	agtcattggga	gtcctgttag	gagtagggat	ggcagagagc	ccccatcagg	tctttaatgt
	aacctggaga	gtcaccaacc	tgatgactgg	gcgtaccgcc	aatgccacct	ccctcctggg

EP 2 913 402 B1

	aactgtacaa	gatgccttcc	caaaattata	ttttgatcta	tgtgatctgg	tcggagagga
	gtgggaccct	tcagaccagg	aaccgtatgt	cgggtatggc	tgcaagtacc	ccgcagggag
	acagcggacc	cggacttttg	acttttacgt	gtgccctggg	cataccgtaa	agtcgggggtg
5	tgggggacca	ggagagggct	actgtggtaa	atgggggtgt	gaaaccaccg	gacaggctta
	ctggaagccc	acatcatcgt	gggaccta	ctcccttaag	cgcgtaaca	ccccctggga
	cacgggatgc	tctaaagtgt	cctgtggccc	ctgctacgac	ctctccaaag	tatccaattc
	cttccaaggg	gctactcgag	ggggcagatg	caaccctcta	gtcctagaat	tcactgatgc
	aggaaaaaag	gctaactggg	acgggccc	atcgtgggga	ctgagactgt	accggacagg
	aacagatcct	attaccatgt	tctccctgac	ccggcaggtc	cttaatgtgg	gaccccgagt
10	ccccataggg	cccaaccag	tattaccgga	ccaaagactc	ccttcctcac	caatagagat
	tgtaccggct	ccacagccac	ctagccccct	caataccagt	tacccccctt	ccactaccag
	tacaccctca	acctcccc	caagtccaag	tgtcccacag	ccacccccag	gaactggaga
	tagactacta	gctctagtca	aaggagccta	tcaggcgctt	aacctcacca	atccccgaca
	gacccaagaa	tgttggtgtg	gcttagtgct	gggacctcct	tattacgaag	gagtagcggt
	cgtgggcact	tataccaatc	attccaccgc	tccggccaac	tgtacggcca	cttcccaaca
15	taagcttacc	ctatctgaag	tgacaggaca	gggcctatgc	atgggggcag	tacctaaaac
	tcaccaggcc	ttatgtaaca	ccacccaaag	cgccggctca	ggatcctact	accttgacgc
	accgcgcgga	acaatgtggg	cttgacgac	tggattgact	ccctgcttgt	ccaccacggt
	gctcaatcta	accacagatt	attgtgtatt	agttgaactc	tggcccagag	taattttacca
	ctcccccgat	tatatgtatg	gtcagcttga	acagcgtacc	aaatataaaa	gagagccagt
20	atcattgacc	ctggcccttc	tactaggagg	attaaccatg	ggagggattg	cagctggaat
	agggacgggg	accactgcct	taattaaaac	ccagcagttt	gagcagcttc	atgcccgtat
	ccagacagac	ctcaacgaag	tcgaaaagtc	aattaccaac	ctagaaaagt	cactgacctc
	gttgtctgaa	gtagtcttac	agaaccgcag	aggcctagat	ttgctattcc	taaaggaggg
	aggctctctgc	gcagccctaa	aagaagaatg	ttgtttttat	gcagaccaca	cggggctagt
	gagagacagc	atggccaaat	taagagaaag	gcttaatcag	agacaaaaac	tatttgagac
25	aggccaagga	tggttcgaag	ggctgtttta	tagatcccc	tggtttacca	ccttaatctc
	caccatcatg	ggacctctaa	tagtactctt	actgatctta	ctctttggac	cttgcatctc
	caatcgattg	gtccaatttg	ttaaagacag	gatctcagtg	gtccaggctc	tggttttgac
	tcagcaatat	caccagctaa	aaccataga	gtacgagcca	tgaacgcgtt	actggccgaa
	gccgcttgga	ataaggccgg	tgtgcgtttg	tctatatgtt	attttccacc	atattgccgt
30	cttttgga	tgtgagggcc	cggaaacctg	gccctgtctt	cttgacgagc	attcctaggg
	gtctttcccc	tctcgccaaa	ggaatgcaag	gtctgttgaa	tgtcgtgaag	gaagcagttc
	ctctggaagc	ttcttgaa	caaacaacgt	ctgtagcgac	cctttgcagg	cagcggaa
	ccccacctgg	cgacaggtgc	ctctgcggcc	aaaagccacg	tgtataagat	acacctgcaa
	aggcggcaca	acccagtg	cacgttgtga	gttgatagt	tgtggaaga	gtcaaattggc
35	tctcctcaag	cgtattcaac	aaggggctga	aggatgcca	gaaggtacct	taagtgtagg
	gatctgatct	ggggcctcgg	tgacatgct	ttacatgtgt	ttagtcgagg	ttaaaaaac
	gtctaggccc	cccgaaccac	ggggacgtgg	ttttcctttg	aaaaacacga	ttataaatgg
	tgaccggcgg	catggcctcc	aagtgggatc	aaaagggcat	ggatatcgct	tacgaggagg
	cctgtctggg	ctacaaggag	ggcggcgtgc	ctatcggcgg	ctgtctgatc	aacaacaagg
	acggcagtg	gctgggcagg	ggccacaaca	tgaggttcca	gaagggtccc	gccacctgc
40	acggcgagat	ctccacctg	gagaactgtg	gcaggctgga	gggcaagggtg	tacaaggaca
	ccacctgtga	caccacctg	tccccttgtg	acatgtgtac	cggcgctatc	atcatgtacg
	gcatccctag	gtgtgtgatc	ggcgagaacg	tgaacttcaa	gtccaagggc	gagaagtacc
	tgcaaaccag	gggccacgag	gtggtggttg	ttgacgatga	gaggtgtaag	aagctgatga
	agcagttcat	cgacgagagg	cctcaggact	ggttcgagga	tatcgtcgag	taagcggccg
45	cagataaaat	aaaagatttt	atttagtctc	cagaaaaagg	ggggaatgaa	agacccacc
	tgtaggtttg	gcaagctagc	ttaagtaacg	ccattttgca	aggcatggaa	aaatacataa
	ctgagaatag	agaagttcag	atcaaggtca	ggaacagatg	gaacagctga	atatgggcca
	aacaggatat	ctgtggtaag	cagttcctgc	cccggctcag	ggccaagaac	agatggaaca
	gctgaatatg	ggccaaacag	gatatctgtg	gtaagcagtt	cctgccccgg	ctcagggcca
	agaacagatg	gtccccagat	gcggtccagc	cctcagcagt	ttctagagaa	ccatcagatg
50	tttccagggt	gccccaaagg	cctgaaatga	ccctgtgcct	tattttgaact	aaccaatcag
	ttcgcttctc	gcttctgttc	gcgcgcttct	gtcccccgag	ctcaataaaa	gagcccacaa
	ccccctacte	ggggcgccag	tctccgatt	gactgagtcg	cccgggtacc	cgtgtatcca
	ataaaccttc	ttgcagttgc	atccgacttg	tggctctcgt	gttccttggtg	agggctcct
	ctgagtgatt	gactaccgt	cagcgggggt	ctttcattac	atgtgagcaa	aaggccagca
55	aaaggccagg	aaccgtaaaa	aggccgcgtt	gctggcgttt	ttccataggc	tccgcccccc
	tgacgagcat	cacaaaaatc	gacgctcaag	tcagaggtgg	cgaaaccgga	caggactata
	aagataccag	gcgtttcccc	ctggaagctc	cctcgtgcgc	tctcctgttc	cgacctgcc

```

gcttaccgga tacctgtccg cctttctccc ttccgggaagc gtggcgcttt ctcatagctc
acgctgtagg tatctcagtt cgggtgtaggt cggtcgctcc aagctgggct gtgtgcacga
acccccggtt cagcccgacc gctgcgcctt atccggtaac tatcgtcttg agtccaaccc
5 ggtaagacac gacttatcgc cactggcagc agccactggt aacaggatta gcagagcgag
gtatgtaggc ggtgctacag agttcttgaa gtgggtggct aactacggct acactagaag
gacagtatth ggtatctgcg ctctgctgaa gccagttacc ttccgaaaaa gagttggtag
ctcttgatcc ggcaaaaaa ccaccgctgg tagcgggtgt ttttttgttt gcaagcagca
gattacgcgc agaaaaaaag gatctcaaga agatcctttg atcttttcta cggggtctga
cgctcagtg aacgaaaact caggttaagg gatttttggtc atgagattat caaaaaggat
10 cttcacctag atccttttaa attaaaaatg aagtttttaa tcaatctaaa gtatatatga
gtaaacttgg tctgacagtt accaatgctt aatcagttag gcacctatct cagcgatctg
tctatttctg tcatccatag ttgcctgact ccccgctgtg tagataacta cgatacggga
gggcttacca tctggcccca gtgctgcaat gataccgcga gaccacgct caccggctcc
agatttatca gcaataaacc agccagccgg aaggggcgag cgcagaagtg gtctgcaac
15 tttatccgcc tccatccagt ctattaattg ttgcccggaa gctagagtaa gtagttcgcc
agttaatagt ttgcgcaacg ttgttgccat tgctgcaggc atcgtggtgt cacgctcgtc
gtttggtatg gcttcattca gctccggttc ccaacgatca aggcgagtta catgatcccc
catgttggtgc aaaaaagcgg ttagctcctt cggctcctcg atcgttggtca gaagtaagtt
ggccgcagtg ttatcactca tggttatggc agcactgcat aattctctta ctgtcatgcc
20 atccgtaaga tgcttttctg tgactgggtga gtactcaacc aagtcattct gagaatagtg
tatgcggcga ccgagttgct cttgcccggc gtcaacacgg gataataacc cgccacatag
cagaacttta aaagtgtca tcattggaaa acgttcttcg gggcgaaaaa tctcaaggat
cttaccgctg ttgagatcca gttcgatgta acccactcgt gcacccaact gatcttcagc
atcttttact ttcaccagcg tttctgggtg agcaaaaaa ggaaggcaaa atgccgcaaa
25 aaagggaata agggcgacac ggaaatgttg aatactcata ctcttccttt ttcaatatta
ttgaagcatt tatcagggtt attgtctcat gagcggatag atatttgaat gtatttagaa
aaataaacia ataggggttc cgcgcacatt tccccgaaa gtgccacctg acgtctaaga
aaccattatt atcatgacat taacctataa aaataggcgt atcacgaggc cctttcgtct
tcaagaattc cat

```

[0018] Disclosed herein are assays for detection of target molecules, such as target nucleic acids comprising viral RNA or DNA, using a plurality of nucleic acid amplification techniques including, for example: NAAT (J.D. Fox, J.Clin. Virol. 40 Suppl. 1 S15-S23, 2007), PCR, RT-PCR, qPCR and RT-qPCR with "touchdown" modifications to improve sensitivity to single copy/assay; RNA transcription based assays (e.g. analogous to the HIV-1 Aptiva assay, see <http://www.fda.gov/BiologicsBloodVaccines/BloodBloodProducts/ApprovedProducts/LicensedProductsBLAs/BloodDonorScreening/InfectiousDisease/ucm149922.htm>); the "branched DNA" system (see, e.g., Anastassopoulou et al., Journal of Virological Methods, 91:67-74, 2001); and further variations in NAAT known to those of skill in the art. Also disclosed are methods using NAAT on nucleic acid samples extracted from histologically fixed samples.

[0019] The assays provided by the disclosure may be used in many embodiments to detect sequence-specific nucleic acids. Disclosed herein are different embodiments of assays using amplification (e.g., PCR) and enzymatic degradation of RNA/DNA heteroduplexes.

[0020] Generally, the disclosure provides a method of identifying MLV-related viral polynucleotides in a subject or sample. The disclosure utilizes a combination of primers and probes having identity to conserved regions of MLV-related viruses. The primers are used to amplify target polynucleotides in the sample and probes are then used to visualize or detect the amplified products. Typically the probe is detectably labeled for detection (e.g., fluorescently labeled, luminescently labeled, enzyme conjugated, radionucleotide labeled and the like). One advantage of the disclosure is that the primer pairs can be used to amplify MLV-related polynucleotides in a sample such as MLV, XMRV or MLV and XMRV polynucleotides. This is advantageous for the detection of XMRV and naturally occurring variants thereof as well as for detecting MLV and naturally occurring variants thereof (including recombinantly engineered MLV vectors).

[0021] For example, a combination of primers and probes identified herein can be used to identify or detect XMRV in a sample, tissue or subject by using primer pairs having homology to MLV (e.g., primer pairs that share at least 95% sequence identity between and XMRV viral sequence and an MLV viral sequence) and a probe sequence that is specific for XMRV (e.g., the probe only hybridizes to an amplified product under highly stringent conditions). Primers that share such homology between XMRV and MLV are identified in Figure 5 and Table 2.

[0022] One utility of this general method is to screen blood and tissue supplies for infection. XMRV has been suggested to be associated with prostate cancer and chronic fatigue syndrome. In another utility, a subject may be screened prior to undergoing treatment with a recombinant retroviral vector. By identifying subject that may have circulating viral polynucleotides in the tissue a risk of recombination between the inherent viral polynucleotide and the therapeutic viral polynucleotide may be managed.

[0023] In another general example, a combination of primers and probes identified herein can be used to identify or detect MLV or MLV-related polynucleotides in a sample, tissue or subject by using primer pairs having homology to MLV (e.g., primer pairs that share at least 95% sequence identity between and XMRV viral sequence and an MLV viral sequence) and a probe sequence that is specific for MLV (e.g., the probe only hybridizes to an amplified product under highly stringent conditions).

[0024] In another embodiment of the general methods described herein, the methods provide useful diagnostics for monitoring patients after delivery of a replication competent MLV-related viral vector. The method can be used to monitor a subject following delivery of the vector on a routine basis (e.g., weekly, monthly, yearly) for as long as a treating physician deems necessary.

[0025] As used herein "MLV-related virus" refers to a retrovirus comprising the general structure of an MLV virus (e.g., LTR-gag-pol-env-LTR) and having at least 60% identity to any of the following sequences set forth in the identified accession numbers: Xenotropic murine leukemia virus isolate LAPC4, complete genome (8,657 bp linear DNA, JF908816.1 GI:336462519); Xenotropic murine leukemia virus isolate VCaP, complete genome (8,657 bp linear DNA, JF908815.1 GI:336462515); Murine leukemia virus N417, complete genome (8,189 bp linear RNA, HQ246218.1 GI:313762331); Moloney murine leukemia virus neuropathogenic variant ts1-92b, complete genome (8,332 bp linear DNA, AF462057.1 GI:18448741); DG-75 Murine leukemia virus, complete genome (8,207 bp linear RNA, AF221065.1 GI:11078528); Rauscher murine leukemia virus, complete genome (8,282 bp linear DNA, NC_001819.1 GI:9629514); Murine leukemia virus SL3-3, complete genome (8,377 bp linear RNA, AF169256.1 GI:5881088); Murine leukemia virus strain SRS 19-6 complete genome (8,256 bp linear DNA, AF019230.1 GI:4071074); Mus dunni endogenous virus complete genome (8,655 bp linear DNA, AF053745.1 GI:3309122); Moloney murine leukemia virus, complete genome (8,332 bp linear RNA, AF033811.1 GI:2801468); Murine type C retrovirus, complete genome (8,135 bp linear DNA, NC_001702.1 GI:9628654); Rauscher murine leukemia virus, complete genome (8,282 bp linear DNA, U94692.1 GI:2228757); Murine leukemia virus isolate NeRV, complete genome (8,273 bp linear RNA, DQ366149.1 GI:86651892); Murine leukemia virus serotype HEMV provirus, complete genome (8,546 bp linear DNA, AY818896.1 GI:55979252); Murine leukemia virus strain BM5eco, complete genome (8,281 bp linear DNA, AY252102.1 GI:30908470); Murine leukemia virus MCF1233, complete genome (8,196 bp linear DNA, U13766.1 GI:535516); MuLV (strain RadLV/VL3(T+L+)) RNA, complete genome (8,394 bp linear RNA, K03363.1 GI:332032); Mink cell focus-forming 247 MuLV env gene, 3' end and LTR (1,164 bp linear RNA, J02249.1 GI:332023); Friend murine leukemia virus, complete genome (8,282 bp linear RNA, M93134.1 GI:331898); Friend murine leukemia virus, complete genome (8,323 bp linear RNA, NC_001362.1 GI:9626096); Friend murine leukemia virus (F-MuLV) complete RNA genome (8,359 bp linear RNA, X02794.1 GI:61544); Gallus gallus MLV-related endogenous retrovirus, complete genome (9,133 bp linear DNA, DQ280312.2 GI:169805278); Friend murine leukemia virus genomic RNA, complete genome, clone:A8 (8,358 bp linear RNA, D88386.1 GI:2351211); PreXMRV-1 provirus, complete genome (8,197 bp linear DNA; NC_007815.2 GI:339276104); Xenotropic MuLV-related virus RKO, complete genome (8,172 bp linear DNA, JF274252.1 GI:338191621); XMRV complete proviral genome, isolate S-162 (8,562 bp linear DNA, FR872816.1 GI:336087897); PreXMRV-2 complete proviral genome (8,193 bp linear DNA, FR871850.1 GI:334849718); Xenotropic MuLV-related virus 22Rv1/CWR-R1 complete proviral genome (8,185 bp linear DNA, FN692043.2 GI:334717372); Xenotropic MuLV-related virus isolate xmlv15, complete genome (8,176 bp linear RNA, HQ154630.1 GI:320091412); PreXMRV-1 complete proviral genome (8,197 bp linear DNA, FR871849.1 GI:334849715); Xenotropic MuLV-related virus VP62, complete genome (8,185 bp linear RNA, DQ399707.1 GI:88765817); Xenotropic MuLV-related virus VP42, complete genome (8,185 bp linear RNA, DQ241302.1 GI:82582299); Xenotropic MuLV-related virus VP35, complete genome (8,185 bp linear RNA, DQ241301.1 GI:82582295); Xenotropic MuLV-related virus VP62, complete genome (8,165 bp linear RNA, EF185282.1 GI:121104176); Plasmid pAMS with hybrid amphotropic/Moloney murine leukemia virus, complete sequence (11,328 bp circular DNA, AF010170.1 GI:2281586); Amphotropic murine leukemia virus strain 1313, complete genome (8,217 bp linear DNA, AF411814.1 GI:28892668); Toca511, recombinant replication competent MLV comprising a polynucleotide encoding cytosine deaminase (see, e.g., SEQ ID NO:19, 20 and 22 of PCT/US2009/058512).

[0026] Any number of different alignment programs can be used to identified regions of identity between any combination of the foregoing MLV-related genomes. Other genomes will be readily identified by using a BLAST algorithm or other similar algorithm to identify sequences having homology/identity to the foregoing sequences.

[0027] In another embodiment, the disclosure relates to a method of detecting XMRV or other MLV related nucleic acids in a sample by using primers and probes that are not specific to XMRV but rather are shared between XMRV and other related strains of MLV. In a related embodiment the MLV/XMRV specific primers are SEQ ID NOs:7 and 8 and the probe is SEQ ID NO:9 (Pol 2 primers and probe; other primers and probes are set forth in Table 1, 2, 3, and 4 below). In a further related embodiment, the MLV/XMRV specific primers are SEQ ID NOs:10 and 11 and the probe is SEQ ID NO:12 (pol1 primer and probe). In a further related embodiment, other sequences can be identified that are common to XMRV and MLV (see the BLAST sequence comparison of two genomes of XMRV and MLV, Figure 5, where perfect sequence homologies of 20 or more bases are underlined/highlighted). Such homologous sequences (or shorter runs of homology down to 15 bases) can be used to select primers and probes. Alternatively, primers and probes can be

chosen using programs that compare sequences and suggest common primers and probes. Such programs are usually designed to look for related genes in different species, and can be used to design QPCR reagents for molecules such as MLV and XMRV with significant but incomplete homology. An example of such a program is Primaclade <http://www.umsl.edu/services/kellogg/primaclade/FAQ.html>.

[0028] An oligonucleotide probe or a primer refers to a nucleic acid molecule of between 8 and 2000 nucleotides in length, or is specified to be about 6 and 1000 nucleotides in length. More particularly, the length of these oligonucleotides can range from about 8, 10, 15, 20, or 30 to 100 nucleotides, but will typically be about 10 to 50 (e.g., 15 to 30 nucleotides). The appropriate length for oligonucleotides in assays of the disclosure under a particular set of conditions may be empirically determined by one of skill in the art. As used herein a primer or probe consisting of at least 95% identity to a reference sequence means that the sequence comprises a sequence that is the same number of oligonucleotides in length, but may differ in nucleotides by 5% from the reference sequence. In addition, a primer or probe consisting of at least 95% identity to a reference sequence and having from 1-10 additional or deleted nucleotides at the 5' and/or 3' end of the oligonucleotide means that the sequence differs by 1 to up to 20 nucleotides in length from the reference sequence and which also differs 5% or less in identity. Accordingly, any primer probe disclosed herein by consist of a reference sequence (e.g., SEQ ID NO:1-21 or sequences set forth in Table 1 and 2); can consist of a sequence that is 95% of greater in identity to a reference sequence (e.g., SEQ ID NO:1-21 or sequences set forth in Table 1 and 2); or can consist of a sequence that is at least 95% identical and has an additional 1-20 (e.g., 1, 2, 3, 4, 5, 6, 7, 8, 9 or 10 nucleotides at the 5' and/or 3' end of the reference sequence (e.g., SEQ ID NO:1-21 or sequence set forth in Table 1 and 2).

[0029] Oligonucleotide primers and probes can be prepared by any suitable method, including direct chemical synthesis and a number of probe systems with derivatized oligonucleotides are available to hybridize to and to detect amplified product, normally by fluorescence change on binding to the probe target. The oligonucleotide primers and probes can contain conventional nucleotides, as well as any of a variety of analogs. For example, the term "nucleotide", as used herein, refers to a compound comprising a nucleotide base linked to the C-1' carbon of a sugar, such as ribose, arabinose, xylose, and pyranose, and sugar analogs thereof. The term nucleotide also encompasses nucleotide analogs. The sugar may be substituted or unsubstituted. Substituted ribose sugars include, but are not limited to, those riboses in which one or more of the carbon atoms, for example the 2'-carbon atom, is substituted with one or more of the same or different Cl, F, --R, --OR, --NR₂ or halogen groups, where each R is independently H, C₁-C₆ alkyl or C₅-C₁₄ aryl. Exemplary riboses include, but are not limited to, 2'-(C₁-C₆)alkoxyribose, 2'-(C₅-C₁₄)aryloxyribose, 2',3'-didehydroribose, 2'-deoxy-3'-haloribose, 2'-deoxy-3'-fluororibose, 2'-deoxy-3'-chlororibose, 2'-deoxy-3'-aminoribose, 2'-deoxy-3'-(C₁-C₆)alkylribose, 2'-deoxy-3'-(C₁-C₆)alkoxyribose and 2'-deoxy-3'-(C₅-C₁₄)aryloxyribose, ribose, 2'-deoxyribose, 2',3'-dideoxyribose, 2'-haloribose, 2'-fluororibose, 2'-chlororibose, and 2'-alkylribose, e.g., 2'-O-methyl, 4'- α -anomeric nucleotides, 1'- α -anomeric nucleotides, 2'-4'- and 3'-4'-linked and other "locked" or "LNA", bicyclic sugar modifications (see, e.g., PCT published application nos. WO 98/22489, WO 98/39352; and WO 99/14226).

[0030] Modifications at the 2'- or 3'-position of ribose include, but are not limited to, hydrogen, hydroxy, methoxy, ethoxy, allyloxy, isopropoxy, butoxy, isobutoxy, methoxyethyl, alkoxy, phenoxy, azido, amino, alkylamino, fluoro, chloro and bromo. Nucleotides include, but are not limited to, the natural D optical isomer, as well as the L optical isomer forms (see, e.g., Garbesi (1993) Nucl. Acids Res. 21:4159-65; Fujimori (1990) J. Amer. Chem. Soc. 112:7435; Urata, (1993) Nucleic Acids Symposium Ser. No. 29:69-70). When the nucleotide base is purine, e.g. A or G, the ribose sugar is attached to the N₉-position of the nucleotide base. When the nucleotide base is pyrimidine, e.g. C, T or U, the pentose sugar is attached to the N₁-position of the nucleotide base, except for pseudouridines, in which the pentose sugar is attached to the C₅ position of the uracil nucleotide base (see, e.g., Kornberg and Baker, (1992) DNA Replication, 2nd Ed., Freeman, San Francisco, Calif.). The 3' end of the probe can be functionalized with a capture or detectable label to assist in detection of a target polynucleotide or of a polymorphism.

[0031] Any of the oligonucleotides or nucleic acids of the disclosure can be labeled by incorporating a detectable label measurable by spectroscopic, photochemical, biochemical, immunochemical, or chemical means. For example, such labels can comprise radioactive substances (e.g., ³²P, ³⁵S, ³H, ¹²⁵I), fluorescent dyes (e.g., 5-bromodesoxyuridin, fluorescein, acetylaminofluorene, digoxigenin), molecular quenchers (e.g. blackhole, molecular beacons), biotin, nanoparticles, and others known to those skilled in the art. Such oligonucleotides are typically labeled at their 3' and/or 5' ends.

[0032] A probe refers to a molecule which can detectably distinguish changes in gene expression or can distinguish between target molecules differing in structure. Detection can be accomplished in a variety of different ways depending on the type of probe used and the type of target molecule. Thus, for example, detection may be based on discrimination of activity levels of the target molecule, but typically is based on detection of specific binding. Examples of such specific binding include nucleic acid probe hybridization. Thus, for example, probes can include nucleic acid hybridization probes (including primers useful for polynucleotide amplification and/or detection). Thus, in one embodiment, the detection of the presence or absence of the at least one target polynucleotide involves contacting a biological sample with a probe or primer pair. Typically an oligonucleotide probe or primer pair, where the probe/primers hybridizes with a form of a target polynucleotide in the biological sample containing a complementary sequence, undergoes hybridization using hybridization selective conditions. Such an oligonucleotide probe can include one or more nucleic acid analogs, labels

or other substituents or moieties so long as the base-pairing function is retained.

[0033] The disclosure provides methods and systems for identifying and quantifying the amount of a given nucleic acid sequence in a given sample, usually down to a single copy per sample. Furthermore, the methods and systems of the disclosure provide sequence specific detection useful for differentiating/identifying related genomic sequences and provide a detectable signal when the correct target sequence is present.

[0034] Methods known in the art can be used to quantitatively measure the amount of a nucleic acid present in a sample. Examples of such methods include quantitative polymerase chain reaction (qPCR), and other NAAT technologies as described above.

[0035] In one embodiment, a method for detecting a specific viral polynucleotide is provided by the disclosure. Such a method can include the use of primers, probes, enzymes, and other reagents for the preparation, detection, and quantitation of a viral polynucleotide (e.g., by PCR, by Northern blot and the like). The primers listed in SEQ ID NOs: 7-12 are particularly suited for use in profiling using RT-PCR based on a viral polynucleotide. Other primers (sense and antisense) and probes are provided in Table 2. One of skill in the art can readily identify in Table 2, the combinations of primer/probes useful for detection of an MLV, MLV-related or XMRV viral sequence. Although the disclosure provides particular primers and probes, those of skill in the art will readily recognize that additional probes and primers can be generated based upon the polynucleotide sequences provided by the disclosure (see, also, for example, Figure 5). Referring to the primers and probes exemplified herein, a series of primers were designed to amplify portions of a murine retroviral (MLV) genome. The primer/probe sets listed in SEQ ID NOs: 1-12 were designed, selected, and tested accordingly (see Examples). Though a number of detection schemes for detecting amplicons are contemplated, as will be discussed in more detail below, one method for detection of polynucleotide amplicons is fluorescence spectroscopy, and therefore labels suited to fluorescence spectroscopy are desirable for detecting polynucleotide. In a related form of detection the amplicon polynucleotide is detected without a hybridization probe but directly with a fluorophore that binds DNA and fluoresces at a wavelength different from that of the free reagent. An example of such a fluorescent label is SYBR Green, though numerous related fluorescent molecules are known including, without limitation, DAPI, Cy3, Cy3.5, Cy5, CyS.5, Cy7, umbelliferone, fluorescein, fluorescein isothiocyanate (FITC), rhodamine, dichlorotriazinylamine fluorescein, dansyl chloride or phycoerythrin.

[0036] In one embodiment of the disclosure, an oligonucleotide primer pair is used to amplify a polynucleotide corresponding to the *pol* gene of a murine retrovirus. Primers comprising SEQ ID NOs:7 and 8 (forward and reverse primers, respectively) are used to amplify the region of the *pol* gene. The amplified region can then be detected using a probe (SEQ ID NO:9), which specifically hybridizes to the amplified polynucleotide. The probe can be labeled with any number of detectable labels as described herein.

[0037] The foregoing primers (SEQ ID NOs:7, 8 and 10, 11) with the appropriate probes (SEQ ID NOs:9 and 12, respectively) can be used to detect the presence of, for example, Murine Leukemia Virus (MLV) and Xenotropic Murine Retrovirus (XMRV). As described elsewhere herein identifying the presence of MLV and/or XMRV is useful for the determination of a therapeutic gene delivery and cancer treatment regimen.

[0038] There is also described an oligonucleotide primer pair is used to amplify a polynucleotide corresponding to the *gag* gene of a XMRV. Primers comprising SEQ ID NOs: 1 and 2 (forward and reverse primers, respectively) are used to amplify the region of the *gag* gene of XMRV. The amplified region can then be detected using a probe (SEQ ID NO:3) which specifically hybridizes to the amplified polynucleotide. The probe can be labeled with any number of detectable labels as described herein. This combination of primers and probes is useful for specifically identifying the presence of an XMRV infection or contamination.

[0039] In another embodiment of the disclosure, an oligonucleotide primer pair is used to amplify a polynucleotide corresponding to the LTR of an MLV vector. Primers comprising SEQ ID NOs: 16 and 17 (forward and reverse primers, respectively) are used to amplify the region of the LTR of an MLV vector (e.g., Toca511). The amplified region can then be detected using a probe (SEQ ID NO:18) which specifically hybridizes to the amplified polynucleotide. The probe can be labeled with any number of detectable labels as described herein. This combination of primers and probes is useful for specifically identifying the presence of a retroviral vector during gene delivery monitoring.

[0040] In another embodiment of the disclosure, an oligonucleotide primer pair is used to amplify a polynucleotide corresponding to a polynucleotide encoding a cytosine deaminase. Primers comprising SEQ ID NOs:19 and 20 (forward and reverse primers, respectively) are used to amplify the region of the a cytosine deaminase delivered using an MLV vector (e.g., Toca511). The amplified region can then be detected using a probe (SEQ ID NO:21) which specifically hybridizes to the amplified polynucleotide. The probe can be labeled with any number of detectable labels as described herein. This combination of primers and probes is useful for specifically identifying the presence of a retroviral vector during gene delivery monitoring.

[0041] The primers (SEQ ID NOs: 1 and 2) and probe (SEQ ID NO:3) can be used to detect the presence of XMRV. As described elsewhere herein identifying the presence of XMRV is useful for the determination of a therapeutic gene delivery, cancer treatment regimen and blood supply screening.

[0042] There is also described an oligonucleotide primer pair is used to amplify a polynucleotide corresponding to the

env gene of a XMRV. Primers comprising SEQ ID NOs:4 and 5 (forward and reverse primers, respectively) are used to amplify the region of the *env* gene of XMRV. The amplified region can then be detected using a probe (SEQ ID NO:6) which specifically hybridizes to the amplified polynucleotide. The probe can be labeled with any number of detectable labels as described herein. This combination of primers and probes is useful for specifically identifying the presence of an XMRV infection or contamination.

[0043] The foregoing primers (SEQ ID NOs:4 and 5) and probe (SEQ ID NO:6) can be used to detect the presence of XMRV. As described elsewhere herein identifying the presence of XMRV is useful for the determination of a therapeutic gene delivery, cancer treatment regimen and screening the blood supply.

[0044] Any of the oligonucleotide primers and probes of the disclosure can be immobilized on a solid support. Solid supports are known to those skilled in the art and include the walls of wells of a reaction tray, test tubes, polystyrene beads, magnetic beads, nitrocellulose strips, membranes, microparticles such as latex particles, glass and the like. The solid support is not critical and can be selected by one skilled in the art. Thus, latex particles, microparticles, magnetic or non-magnetic beads, membranes, plastic tubes, walls of microtiter wells, glass or silicon chips and the like are all suitable examples. Suitable methods for immobilizing oligonucleotides on a solid phase include ionic, hydrophobic, covalent interactions and the like. The solid support can be chosen for its intrinsic ability to attract and immobilize the capture reagent. The oligonucleotide probes or primers of the disclosure can be attached to or immobilized on a solid support individually or in groups of about 2-10,000 distinct oligonucleotides of the disclosure to a single solid support.

[0045] A substrate comprising a plurality of oligonucleotide primers or probes of the disclosure may be used either for detecting or amplifying targeted sequences. The oligonucleotide probes and primers of the disclosure can be attached in contiguous regions or at random locations on the solid support. Alternatively the oligonucleotides of the disclosure may be attached in an ordered array wherein each oligonucleotide is attached to a distinct region of the solid support which does not overlap with the attachment site of any other oligonucleotide. Typically, such oligonucleotide arrays are "addressable" such that distinct locations are recorded and can be accessed as part of an assay procedure. The knowledge of the location of oligonucleotides on an array make "addressable" arrays useful in hybridization assays. For example, the oligonucleotide probes can be used in an oligonucleotide chip such as those marketed by Affymetrix and described in U.S. Pat. No. 5,143,854; PCT publications WO 90/15070 and 92/10092. These arrays can be produced using mechanical synthesis methods or light directed synthesis methods which incorporate a combination of photolithographic methods and solid phase oligonucleotide synthesis.

[0046] The immobilization of arrays of oligonucleotides on solid supports has been rendered possible by the development of a technology generally referred to as "Very Large Scale Immobilized Polymer Synthesis" in which probes are immobilized in a high density array on a solid surface of a chip (see, e.g., U.S. Patent Nos. 5,143,854; and 5,412,087 and in PCT Publications WO 90/15070, WO 92/10092 and WO 95/11995), which describe methods for forming oligonucleotide arrays through techniques such as light-directed synthesis techniques.

[0047] In another embodiment, an array of oligonucleotides complementary to subsequences of the target gene (for example yeast cytosine deaminase (CD) or a version of CD optimized for expression in human cells) is used to determine the identity of the target, measure its amount and the like.

[0048] Hybridization techniques can also be used to identify the viral polynucleotides in a subject or sample and thereby determine or predict cross reactivity, chances of recombination or a treatment regimen using a gene delivery vector comprising a recombinant MLV vector. The hybridization reactions may be carried out in a solid support (e.g., membrane or chip) format, in which, for example, a probe (e.g., SEQ ID NO:3, 6 and/or 9) are immobilized on nitrocellulose or nylon membranes and probed with amplified preparations of nucleic acids obtained, for example, from PCR using primers comprising SEQ ID NO:1, 2, 4, 5, 7, 8, 10, 11, 13 and/or 14 of the disclosure. Any of the known hybridization formats may be used, including Southern blots, slot blots, "reverse" dot blots, solution hybridization, solid support based sandwich hybridization, bead-based, silicon chip-based and microtiter well-based hybridization formats.

[0049] Hybridization of an oligonucleotide probe to a target polynucleotide may be performed with both entities in solution, or such hybridization may be performed when either the oligonucleotide or the target polynucleotide is covalently or noncovalently affixed to a solid support. Attachment to a solid support may be mediated, for example, by antibody-antigen interactions, poly-L-Lysine, streptavidin or avidin-biotin, salt bridges, hydrophobic interactions, chemical linkages, UV cross-linking baking, etc. Oligonucleotides may be synthesized directly on the solid support or attached to the solid support subsequent to synthesis. The solid support may be treated, coated or derivatized to facilitate the immobilization of the specific oligonucleotide.

[0050] Hybridization assays based on oligonucleotide arrays rely on the differences in hybridization stability of short oligonucleotides to perfectly matched and mismatched target variants. Each DNA chip can contain thousands to millions of individual synthetic DNA probes arranged in a grid-like pattern and miniaturized to the size of a dime or smaller. Such a chip may comprise oligonucleotides representative of both a wild-type and variant sequences.

[0051] Oligonucleotides of the disclosure can be designed to specifically hybridize to a target region of a polynucleotide. As used herein, specific hybridization means the oligonucleotide forms an anti-parallel double-stranded structure with the target region under certain hybridizing conditions, while failing to form such a structure when incubated with a different

target polynucleotide or another region in the polynucleotide or with a polynucleotide lacking the desired locus under the same hybridizing conditions. Typically, the oligonucleotide specifically hybridizes to the target region under conventional high stringency conditions.

[0052] A nucleic acid molecule such as an oligonucleotide or polynucleotide is said to be a "perfect" or "complete" complement of another nucleic acid molecule if every nucleotide of one of the molecules is complementary to the nucleotide at the corresponding position of the other molecule. A nucleic acid molecule is "substantially complementary" to another molecule if it hybridizes to that molecule with sufficient stability to remain in a duplex form under conventional low-stringency conditions. Conventional hybridization conditions are described, for example, in Sambrook et al., *Molecular Cloning, A Laboratory Manual*, 2nd ed., Cold Spring Harbor Press, Cold Spring Harbor, N. Y. (1989), and in Haymes et al., *Nucleic Acid Hybridization, A Practical Approach*, IRL Press, Washington, D. C. (1985). While perfectly complementary oligonucleotides are used in most assays for detecting target polynucleotides or polymorphisms, departures from complete complementarity are contemplated where such departures do not prevent the molecule from specifically hybridizing to the target region. For example, an oligonucleotide primer may have a non-complementary fragment at its 5' or 3' end, with the remainder of the primer being complementary to the target region. Those of skill in the art are familiar with parameters that affect hybridization; such as temperature, probe or primer length and composition, buffer composition and salt concentration and can readily adjust these parameters to achieve specific hybridization of a nucleic acid to a target sequence.

[0053] A variety of hybridization conditions may be used in the disclosure, including high, moderate and low stringency conditions; see for example Maniatis et al., *Molecular Cloning: A Laboratory Manual*, 2d Edition, 1989, and Short Protocols in *Molecular Biology*, ed. Ausubel, *et al.*. Stringent conditions are sequence-dependent and will be different in different circumstances. Longer sequences hybridize specifically at higher temperatures. An extensive guide to the hybridization of nucleic acids is found in Tijssen, *Techniques in Biochemistry and Molecular Biology--Hybridization with Nucleic Acid Probes*, "Overview of principles of hybridization and the strategy of nucleic acid assays" (1993). Generally, stringent conditions are selected to be about 5-10°C lower than the thermal melting point (T_m) for the specific sequence at a defined ionic strength and pH. The T_m is the temperature (under defined ionic strength, pH and nucleic acid concentration) at which 50% of the probes complementary to the target hybridize to the target sequence at equilibrium (as the target sequences are present in excess, at T_m , 50% of the probes are occupied at equilibrium). Stringent conditions will be those in which the salt concentration is less than about 1.0 M sodium ion, typically about 0.01 to 1.0 M sodium ion concentration (or other salts) at pH 7.0 to 8.3 and the temperature is at least about 30° C for short probes (*e.g.*, 10 to 50 nucleotides) and at least about 60° C for long probes (*e.g.*, greater than 50 nucleotides). Stringent conditions may also be achieved with the addition of helix destabilizing agents such as formamide. The hybridization conditions may also vary when a non-ionic backbone, *i.e.*, PNA is used, as is known in the art. In addition, cross-linking agents may be added after target binding to cross-link, *i.e.*, covalently attach, the two strands of the hybridization complex.

[0054] Methods and compositions of the disclosure are useful for diagnosing or determining the presence of contamination or infection in a sample or subject, respectively. Such tests can be performed using DNA or RNA samples collected from blood, cells, biopsies, tissue scrapings, tissue culture, or other cellular materials. As will be appreciated by those in the art, target polynucleotides can be obtained from samples including, but not limited to, bodily fluids (*e.g.*, blood, urine, serum, lymph, saliva, anal and vaginal secretions, perspiration and semen) of virtually any organism, with mammalian samples common to the methods of the disclosure and human samples being typical. The sample may comprise individual cells, including primary cells (including bacteria) and cell lines including, but not limited to, tumor cells of all types (particularly melanoma, myeloid leukemia, carcinomas of the lung, breast, ovaries, colon, kidney, prostate, pancreas and testes); cardiomyocytes; endothelial cells; epithelial cells; lymphocytes (T-cell and B cell); mast cells; eosinophils; vascular intimal cells; hepatocytes; leukocytes including mononuclear leukocytes; stem cells such as haemopoietic, neural, skin, lung, kidney, liver and myocyte stem cells; osteoclasts; chondrocytes and other connective tissue cells; keratinocytes; melanocytes; liver cells; kidney cells; and adipocytes. Suitable cells also include known research cells, including, but not limited to, Jurkat T cells, NIH3T3 cells, CHO, Cos, 923, HeLa, SiHa, WI-38, Weri-1, MG-63, and the like (see the ATCC cell line catalog).

[0055] Other methods to amplify and identify viral infection or contamination by MLV or XMRV will be recognized in the art and can be utilized in combination with the primers and probes identified herein. For example, one of skill in the art will recognize that Branched DNA, Hybrid Capture Assays, PCR (including RT, nested, multiplex, Real Time), Nucleic acid sequence-based amplification, transcription mediated amplification, strand displacement amplification, Ligase Chain Reaction, Cleavase-invader technology and cycling probe technology can be used with the oligonucleotides of the disclosure.

[0056] A target polynucleotide (*e.g.*, a virus polynucleotide or gene) may be amplified using any oligonucleotide-directed amplification method including, but not limited to, polymerase chain reaction (PCR) (U.S. Pat. No. 4,965,188), ligase chain reaction (LCR) (Barany et al., *Proc. Natl. Acad. Sci. USA* 88:189-93 (1991); WO 90/01069), and oligonucleotide ligation assay (OLA) (Landegren et al., *Science* 241:1077-80 (1988)). Other known nucleic acid amplification procedures may be used to amplify the target region (s) including transcription-based amplification systems (U.S. Pat.

No. 5,130,238; European Patent No. EP 329,822; U.S. Pat. No. 5,169,766; WO 89/06700) and isothermal methods (Walker et al., Proc. Natl. Acad. Sci. USA 89:392-6 (1992)).

[0057] Ligase Chain Reaction (LCR) techniques can be used and are particularly useful for detection of single or multiple (e.g., 1, 2, 3, 4, or 5) nucleotide differences between similar polynucleotides. LCR occurs only when the oligonucleotides are correctly base-paired. The Ligase Chain Reaction (LCR), which utilizes the thermostable Taq ligase for ligation amplification, is useful for interrogating loci of a gene. LCR differs from PCR because it amplifies the probe molecule rather than producing amplicon through polymerization of nucleotides. Two probes are used per each DNA strand and are ligated together to form a single probe. LCR uses both a DNA polymerase enzyme and a DNA ligase enzyme to drive the reaction. Like PCR, LCR requires a thermal cycler to drive the reaction and each cycle results in a doubling of the target nucleic acid molecule. LCR can have greater specificity than PCR. The elevated reaction temperatures permit the ligation reaction to be conducted with high stringency. Where a mismatch occurs, ligation cannot be accomplished. For example, a probe based upon a target polynucleotide is synthesized in two fragments and annealed to the template with possible difference at the boundary of the two primer fragments. A ligase ligates the two primers if they match exactly to the template sequence.

[0058] In one embodiment, the two hybridization probes are designed each with a target specific portion. The first hybridization probe is designed to be substantially complementary to a first target domain of a target polynucleotide (e.g., a polynucleotide fragment) and the second hybridization probe is substantially complementary to a second target domain of a target polynucleotide (e.g., a polynucleotide fragment). In general, each target specific sequence of a hybridization probe is at least about 5 nucleotides long, with sequences of about 15 to 30 being typical and 20 being especially common. In one embodiment, the first and second target domains are directly adjacent, e.g., they have no intervening nucleotides. In this embodiment, at least a first hybridization probe is hybridized to the first target domain and a second hybridization probe is hybridized to the second target domain. If perfect complementarity exists at the junction, a ligation structure is formed such that the two probes can be ligated together to form a ligated probe. If this complementarity does not exist (due to mismatch), no ligation structure is formed and the probes are not ligated together to an appreciable degree. This may be done using heat cycling, to allow the ligated probe to be denatured off the target polynucleotide such that it may serve as a template for further reactions. The method may also be done using three hybridization probes or hybridization probes that are separated by one or more nucleotides, if dNTPs and a polymerase are added (this is sometimes referred to as "Genetic Bit" analysis).

[0059] Quantitative PCR and digital PCR can be used to measure the level of a polynucleotide in a sample. Digital Polymerase Chain Reaction (digital PCR, dPCR or dePCR) can be used to directly quantify and clonally amplify nucleic acids including DNA, cDNA or RNA. Digital PCR amplifies nucleic acids by temperature cycling of a nucleic acid molecule with a DNA polymerase. The reaction is typically carried out in the dispersed phase of an emulsion capturing each individual nucleic acid molecule present in a sample within many separate chambers or regions prior to PCR amplification. A count of chambers containing detectable levels of PCR end-product is a direct measure of the absolute nucleic acids quantity.

[0060] Quantitative polymerase chain reaction (qPCR) is a modification of the polymerase chain reaction and real-time quantitative PCR are useful for measuring the amount of DNA after each cycle of PCR by use of fluorescent markers or other detectable labels. Quantitative PCR methods use the addition of a competitor RNA (for reverse-transcriptase PCR) or DNA in serial dilutions or co-amplification of an internal control to ensure that the amplification is stopped while in the exponential growth phase.

[0061] Modifications of PCR and PCR techniques are routine in the art and there are commercially available kits useful for PCR amplification.

[0062] A probe or primer of the disclosure can be associated with a detectable label. A signaling component can include any label that can be detected optically, electronically, radioactively and the like. A nucleic acid analog may serve as the signaling component. By "label" or "detectable label" is meant a moiety that allows detection. In one embodiment, the detection label is a primary label. A primary label is one that can be directly detected, such as a fluorophore. In general, labels fall into three classes: a) isotopic labels, which may be radioactive or heavy isotopes; b) magnetic, electrical, thermal labels; and c) colored or luminescent dyes. Common labels include chromophores or phosphors but are typically fluorescent dyes. Suitable dyes for use in the disclosure include, but are not limited to; fluorescent lanthanide complexes, including those of Europium and Terbium, fluorescein, rhodamine, tetramethylrhodamine, eosin, erythrosin, coumarin, methyl-coumarins, quantum dots (also referred to as "nanocrystals"), pyrene, Malachite green, stilbene, Lucifer Yellow, Cascade Blue™, Texas Red, Cy dyes (Cy3, Cy5, and the like), Alexa dyes, phycoerythrin, bodipy, and others described in the 6th Edition of the Molecular Probes Handbook by Richard P. Haugland.

[0063] Such a detectable label may be a radioactive label or may be a luminescent, fluorescent or enzyme label. Indirect detection processes typically comprise probes covalently labeled with a hapten or ligand such as digoxigenin (DIG) or biotin. In one embodiment, following the hybridization step, the target-probe duplex is detected by an antibody-streptavidin-enzyme complex. Enzymes commonly used in DNA diagnostics are horseradish peroxidase and alkaline phosphatase. Direct detection methods include the use of fluorophore-labeled oligonucleotides, lanthanide chelate-labeled

oligonucleotides or oligonucleotide-enzyme conjugates. Examples of fluorophor labels are fluorescein, rhodamine and phthalocyanine dyes.

[0064] It will be understood that embodiments of the invention include probes having fluorescent dye molecules, fluorescent compounds, or other fluorescent moieties. A dye molecule may fluoresce, or be induced to fluoresce upon excitation by application of suitable excitation energy (e.g., electromagnetic energy of suitable wavelength), and may also absorb electromagnetic energy ("quench") emitted by another dye molecule or fluorescent moiety. Any suitable fluorescent dye molecule, compound or moiety may be used in the practice of the invention. For example, suitable fluorescent dyes, compounds, and other fluorescent moieties include fluorescein, 6-carboxyfluorescein (6-FAM), 2',4',1,4,-tetrachlorofluorescein (TET), 2',4',5',7',1,4-hexachlorofluorescein (HEX), 2',7'-dimethoxy-4',5'-dichloro-6-carboxyrhodamine (JOE), 2'-chloro-5'-fluoro-7',8'-fused phenyl-1,4-dichloro-6-carboxyfluorescein (NED) and 2'-chloro-7'-phenyl-1,4-dichloro-6-carboxyfluorescein (VIC), cyanine dyes (e.g., Cy.sup.3, Cy.sup.5, Cy.sup.9, nitrothiazole blue (NTB)), Cys3, FAM.TM., tetramethyl-6-carboxyrhodamine (TAMRA), tetrapropano-6-carboxyrhodamine (ROX), dipyrromethene boron fluoride (Bodipy), dichloro-fluorescein, dichloro-rhodamine, fluorescein thiosemicarbazide (FTC), sulforhodamine 101 acid chloride (Texas Red), phycoerythrin, rhodamine, carboxytetramethylrhodamine, 4,6-diamidino-2-phenylindole (DAPI), an indopyras dye, pyrenyloxytrisulfonic acid (Cascade Blue), 514 carboxylic acid (Oregon Green), eosin, erythrosin, pyridyloxazole, benzoxadiazole, aminonaphthalene, pyrene, maleimide, a coumarin, 4-fluoro-7-nitrobenofurazan (NBD), 4-amino-N-[3-(vinylsulfonyl)-phenyl]naphthalimide-3,6-disulfonate (Lucifer Yellow), DABCYL, DABSYL, anthraquinone, malachite green, nitrothiazole, and nitroimidazole compounds, propidium iodide, porphyrins, lanthanide cryptates, lanthanide chelates, derivatives and analogs thereof (e.g., 5-carboxy isomers of fluorescein dyes), and other fluorescent dyes and fluorescent molecules and compounds.

[0065] An oligonucleotide according to the methods of the invention may be labeled at the 5' end or the 3' end of at least one subunit of the probe. In embodiments, oligonucleotides may be labeled at both the 5' end and the 3' end. Alternatively, at least one subunit of the probe may be labeled internally, having at least one, and, in embodiments, more than one, internal label. In embodiments, an oligonucleotide may be labeled at an end and may be labeled internally. The oligonucleotides themselves are synthesized using techniques that are also well known in the art. Methods for preparing oligonucleotides of specific sequence are known in the art, and include, for example, cloning and restriction digest analysis of appropriate sequences and direct chemical synthesis, including, for example, the phosphotriester method described by Narang et al., 1979, *Methods in Enzymology*, 68:190, the phosphodiester method disclosed by Brown et al., 1979, *Methods in Enzymology*, 68:109, the diethylphosphoramidate method disclosed in Beaucage et al., 1981, *Tetrahedron Letters*, 22:1859, and the solid support method disclosed in U.S. Pat. No. 4,458,066, or by other chemical methods using a commercial automated oligonucleotide synthesizer. Modified linkages also may be included, for example phosphorothioates.

[0066] Examples of detection modes contemplated for the disclosed methods include, but are not limited to, spectroscopic techniques, such as fluorescence and UV-Vis spectroscopy, scintillation counting, and mass spectroscopy. Complementary to these modes of detection, examples of labels for the purpose of detection and quantitation used in these methods include, but are not limited to, chromophoric labels, scintillation labels, and mass labels. The expression levels of polynucleotides and polypeptides measured using these methods may be normalized to a control established for the purpose of the targeted determination.

[0067] Label detection will be based upon the type of label used in the particular assay. Such detection methods are known in the art. For example, radioisotope detection can be performed by autoradiography, scintillation counting or phosphor imaging. For hapten or biotin labels, detection is with an antibody or streptavidin bound to a reporter enzyme such as horseradish peroxidase or alkaline phosphatase, which is then detected by enzymatic means. For fluorophor or lanthanide-chelate labels, fluorescent signals may be measured with spectrofluorimeters with or without time-resolved mode or using automated microtitre plate readers. With enzyme labels, detection is by color or dye deposition (p-nitrophenyl phosphate or 5-bromo-4-chloro-3-indolyl phosphate/nitroblue tetrazolium for alkaline phosphatase and 3,3'-diaminobenzidine-NiCl₂ for horseradish peroxidase), fluorescence (e.g., 4-methyl umbelliferyl phosphate for alkaline phosphatase) or chemiluminescence (the alkaline phosphatase dioxetane substrates LumiPhos 530 from Lumigen Inc., Detroit Mich, or AMPPD and CSPD from Tropix, Inc.). Chemiluminescent detection may be carried out with X-ray or polaroid film or by using single photon counting luminometers.

[0068] The methods, compositions, systems and devices disclosed herein find use in the identification and quantization of a target DNA or RNA polynucleotide in a sample, such as in a pool of sequences including one or more target sequences, which may be unrelated polynucleotides. Quantization of specific nucleic acid samples may be achieved by comparing the total signal (fluorescent or otherwise) obtained during the assay with a standard curve of known polynucleotide target concentrations. Specific examples of applications include the detection of pathogenic viruses through the detection of their biomolecules such as DNA or RNA which are indicative of the presence of said targets. Assays disclosed herein may be used to detect and identify the presence of specific DNA sequences and may be used in assays for diagnosis of many types of infection and disease.

[0069] These assays are suitable for use on cell lysates, and contaminated samples as well. Since many clinical

samples are rich with contaminants, it is advantageous that the described assays herein work under these conditions. Although many methods are currently available for DNA extraction and purification from tissues, assays such as those disclosed herein (QIAGEN kits, etc.), which are proficient in analyzing and working with contaminated samples, are very valuable and increases the robustness of the assay. For clinical sample use with the disclosed assays, sample preparation kits may be used. For example, samples suspected of containing pathogenic DNA may be used. Exemplary kits and protocols that can be used include the QIAamp MinElute Virus Spin Kit provided by Qiagen. This kit allows DNA isolation from clinical samples in roughly 1 hour. Other methods for sample preparation are available from suppliers such as Promega.

[0070] Polynucleotides may be prepared from samples using known techniques. For example, the sample may be treated to lyse a cell comprising the target polynucleotide, using known lysis buffers, sonication techniques, electroporation, and the like. Many methods for cell lysis are common knowledge for those trained in the art.

[0071] The following Examples are provided to illustrate and do not limit the invention.

EXAMPLES

[0072] The following abbreviations and definitions will assist in understanding aspect of the disclosure and the assays performed.

[0073] Ct (Cycle Threshold): Cycle number (in qPCR) at which the fluorescence generated within a reaction well exceeds the defined threshold. The threshold is arbitrarily defined by the qPCR instrument manufacturer to reflect the point during the reaction at which a sufficient number of amplicons have accumulated.

gDNA (genomic DNA): Deoxyribonucleic acid that has been purified from tissue and/or cultured cells.

[0074] Percentage Coefficient of Variation (% CV): The coefficient of variation (CV) is a normalized measure of dispersion of a probability distribution. It is defined as the ratio of the standard deviation σ to the mean μ :

$$C_v = \frac{\sigma}{\mu}$$

[0075] Slope: The slope or gradient of a line describes its steepness, incline, or grade. An acceptable slope of the linear regression equation for the qPCR should be within the range of -3.00 to -3.7.

[0076] R-squared (R^2) Value (also known as the Pearson Correlation Coefficient): The correlation of the line, R^2 , is a measure of how well the data fits the model and how well the data fits on a straight line. It is influenced by pipetting accuracy and by the range of the assay. An R^2 of ≥ 0.94 is acceptable.

[0077] qPCR Percentage Efficiency (% Efficiency): Amplification efficiency, E, is the efficiency of amplification at varying template concentrations and is calculated from the slope of the standard curve using the following formula:

$$E = 10^{(-1/\text{slope})}$$

[0078] The % efficiency is the percent of template that was amplified in each cycle and is calculated using the following formula:

$$\% \text{ Efficiency} = (E - 1) \times 100\%$$

The % efficiency should be between 85% and 115%.

LOD: Limit of detection

LLOQ: Lower limit of quantitation.

ND: Non-detected.

NA: Not applicable.

NTC (Non-template control): A series of reaction wells in a qPCR experiment that contains all the reagents necessary for amplification with elution buffer or water substituted for sample DNA.

qPCR: Real-time quantitative polymerase chain reaction.

Example 1: Design of primer/probe sets

[0079] Two primer/probe sets were designed for XMRV specific qPCR:

1. XMRV gag

XMRV 628F (5'-ACTACCCCTCTGAGTCTAACC-3') (SEQ ID NO:1)
 XMRV764R (5'-GGCCATCCTACATTGAAAGTTG-3') (SEQ ID NO:2)
 XMRV gag probe (5'-**FAM**-CGCATTGCATCCAACCACTGTG-3'-BHQ) (SEQ ID NO:3)

Amplification curves are shown in Figure 6.

2. XMRV env

XMRV 6252F (5'-TTTGATTCTCAGTGGGCTC-3') (SEQ ID NO:4)
 XMRV6391R (5'-CGATACAGTCTTAGTCCCCATG-3') (SEQ ID NO:5)
 XMRV env probe (5'-**HEX**-CCCTTTTACCCGCGTCAGTGAATTCT-3'-BHQ) (SEQ ID NO:6)

[0080] Two primer probe sets were designed for detection of all MLV related retroviruses. Amplification curves are shown in Figure 7.

3. MLV Pol 1

pol-F (5'-AACAAGCGGGTGAAGACATC-3') (SEQ ID NO:7)
 pol-R (5'-CAAAGGCGAAGAGAGGCTGAC-3') (SEQ ID NO:8)
 pol probe (5'-**HEX**-CCCACCGTGCCCAACCCCTACAACC-3'-TAMRA) (SEQ ID NO:9)

4. MLV Pol 2

5' Pol2 Primer(CAAGGGGCTACTGGAGGAAAG) (SEQ ID NO:10)
 3' Pol2 Primer:(CTTTCCTCCATGTACCAGACTG) (SEQ ID NO:11)
 Pol2 Probe:(5HEX/TATCGCTGGACCACGGATCGCAA/3BHQ_1) (SEQ ID NO:12)

[0081] Two primer probe sets were designed for detection of amphotropic MLV virus:

5. MLV Env2

5'Env2 primer: 5'-ACCCTCAACCGCCCCTACAAGT-3'(SEQ ID NO: 13)
 3'Env2 primer: 5'-GTTAAGCGCCTGATAGGCTC-3' (SEQ ID NO: 14)
 Env2 probe: 5'-/FAM/CCCCAAATGAAAGACCCCCGCTGACG/BHQ/-3'(SEQ ID NO: 15)

6. MLV LTR:

5' Primer = MLV-U3-B: AGC CCA CAA CCC CTC ACT C (SEQ ID NO:16)
 3' Primer = 3-MLV-Psi: TCT CCC GAT CCC GGA CGA (SEQ ID NO:17)
 FAM Probe = MLV-U5-Psi: FAM-CCCCAAATGAAAGACCCCCGCTGACG 3BHQ_1 (SEQ ID NO:18)

[0082] One primer probe set was designed for detection of a cytosine deaminase (CD) gene.

7. CD:

5' yCD2 Primer: (ATC ATC ATG TAC GGC ATC CCT AG) (SEQ ID NO:19)
 3' yCD2 Primer: (TGA ACT GCT TCA TCA GCT TCT TAC) (SEQ ID NO:20)
 yCD2 Probe: (5FAM/TCA TCG TCA ACA ACC ACC ACC TCG T/3BHQ_1) (SEQ ID NO:21)

[0083] Oligonucleotides for primer probe sets were ordered from IDT (Integrated DNA Technologies, Inc., San Diego, California).

Example 2: Preparation of genomic DNA from blood and other tissues from mammals including human and canines for PCR testing.

[0084] The XMRV (xenotropic murine leukemia virus-related virus) qPCR assay is performed to quantify DNA. Total DNA extraction from the specimens samples is generated by standard means such as the use of commercially available kits (QIAGEN DNA blood mini kit, QIAGEN DNA Tissue kit, Promega DNA Tissue Kit, Promega DNA Cell Kit). A quantitation curve is established with 8 non-zero samples comprising of serial dilutions of defined copy number of reference plasmid to generate a Ct value versus copy number correlation. Linear regression analysis generates an equation which is used to calculate the copy number in the sample. Quantitative curves generation are shown for XMRV gag (Figure 8), XMRV env (Figure 9), XMRV pol2 (Figure 10).

Example 3: Preparation of plasma from humans and dogs for RT-PCR testing.

[0085] Blood was collected in blood collection tubes, and serum or plasma prepared from the whole blood by conventional means. The XMRV (xenotropic murine leukemia virus-related virus) RT-PCR assay is performed to quantify RNA from biological samples, such as whole blood and plasma, without the need for RNA extraction. The assay employs a two-step amplification process with the initial step consisting of the distribution of 2 μ L of experimental sample directly into a cDNA reaction mix. Following completion of the reverse transcriptase (RT) cDNA synthesis, a 2 μ L aliquot is removed, transferred into a qPCR reaction mix and a qPCR protocol is performed. A quantitation curve is established with 6 non-zero samples comprising of serial dilutions of defined copy number of reference vector to generate a Ct value versus copy number correlation. Linear regression analysis generates an equation which is used to calculate the copy number in the sample. Quantitative curves generation are shown in Figure 15.

Example 4: Standardization and Validation of QPCR DNA assays.

[0086] A series of experiments were performed as outlined below:

- 1) To optimize the cycling parameters of the quantitative PCR (qPCR) protocol for XMRV detection including primer and probe concentrations and annealing temperature.
- 2) To assess detection sensitivity in spiked human whole blood genomic DNA (gDNA) targeted with the XMRV env, XMRV gag and XMRV pol2 primer/probe sets using qPCR.
- 3) To assess the use of an additional set of three pre-cycling steps (defined as a stage) in the qPCR protocol with respect to detection sensitivity.
- 4) To assess for variance in XMRV detection sensitivity from independent sources of human whole blood.
- 5) To assess for recovery of 22Rv1 XMRV positive control spiked into human whole blood gDNA.

Assay Design

a. Optimization of cycling parameters for the qPCR protocol.

[0087] A matrix of primers and probe were made up in various concentration combinations and were used to target the appropriate XMRV plasmid (pUC57 XMRV gag, pET28b XMRV env or pAZ3-emd pol2) containing the gene of interest. The choice of optimal primer concentrations were made based on comparisons of Ct value, standard deviation and relative fluorescence units (RFU). SYBR Green was used for the primer concentration optimization qPCR assay and TaqMan was used for the probe concentration optimization assay. Annealing temperature optimization was carried out by performing a qPCR annealing temperature gradient ranging from 50°C to 65°C. Plasmids specific for the gene of interest were targeted with the appropriate XMRV primer sets.

b. Detection sensitivity in human whole blood gDNA spiked with plasmid DNA and targeted with XMRV env, XMRV gag or XMRV pol2 primer/probe sets using qPCR.

[0088] Genomic DNA extracted from human whole blood was spiked with known copy numbers of plasmid DNA containing the gene of interest. Serial log dilutions of the spiked gDNA were made and qPCR was performed. The samples were targeted with XMRV env (Figure 9), XMRV gag (Figure 8) and XMRV pol2 (Figure 10) primer/probe sets in single qPCR reactions.

c. Detection sensitivity in human whole blood gDNA spiked with plasmid DNA and targeted with XMRV env, XMRV gag or XMRV pol2 primer/probe sets using a one-stage qPCR protocol.

[0089] Human whole blood gDNA was spiked with known copy numbers of plasmid DNA containing the gene of interest. Serial log dilutions of the spiked gDNA were made and a modified version of the qPCR protocol was performed by adding a set of three pre-cycling steps (defined as a one-stage qPCR protocol) to the current qPCR protocol. The samples were targeted with XMRV env (Figure 9), XMRV gag (Figure 8) and XMRV pol2 (Figure 10) primer/probe sets in single qPCR reactions.

d. Assessment of XMRV detection sensitivity from human whole blood sourced from healthy donors and spiked with plasmid DNA.

[0090] Genomic DNA from whole blood from healthy donors were used to spike in known copy numbers of plasmid DNA. Serial dilutions of the gDNA were made to generate 1E3, 1E2, 1E1 and 1E0 copies per reaction. A 0-stage and a 1-stage qPCR protocol were performed. The samples were targeted with XMRV env (Figure 12), XMRV gag (Figure 11) or XMRV pol2 (Figure 13) primer/probe sets in single qPCR reactions.

e. 22rev1 Positive control recovery assessment spiked into human whole blood gDNA.

[0091] 22Rv1 gDNA (positive for XMRV, E.C. Knouf et al. J.Virol 83:78353-7356 2009) was spiked into purified human whole blood gDNA (pre and post gDNA extraction) at increasing log dilutions (one human whole blood sample control and one TE sample were spiked pre-extraction with 500 ng of 22Rv1 gDNA to yield a final concentration of ~2.5 ng/μL). A 0-stage and a 1-stage qPCR protocol were performed with primers targeting the XMRV gag, XMRV env and XMRV pol sequences.

[0092] Optimizations (primer concentration and temperature) for Pol primer set were carried out.

[0093] Both XMRV gag and XMRV env primer sets are XMRV specific whereas the Pol primer sets detect both MLV and XMRV.

[0094] A qPCR protocol used for all 4 primer sets: BioRad Supermix65.prc1

Step 1:	95°C 5 min
Step 2:	95°C 15 sec
Step 3:	65°C 30 sec [repeat step 2-3 44X more times]

[0095] Figures 6 and 7 show results obtained by the methods and compositions disclosed above.

[0096] TaqMan® Gold RT-PCR Kit and TaqMan® PCR universal master mix are obtained from PE Biosystems. RNAeasy® mini kit and QIAamp® viral RNA mini kit are obtained from Qiagen. Various cell culture materials and biological samples to be tested are obtained from vendors or subjects.

[0097] MLV recombinant isolates comprise the sequences set forth in International Application No. PCT/US09/58512 and published on April 1, 2010 as publication no. WO 2010/036986.

[0098] Two primer/probe sets for the detection of XMRV were designed as set forth above. One forward primer (FP), one reverse primer (RP), and one probe were used for the detection of XMRV gag and XMRV env. A third set of primer/probe was used for the detection of XMRV and MLV using the primers above that amplify the pol region of XMRV and MLV.

[0099] The qPCR reaction mixture contains 900 nM primers (both forward and reverse) and 200 nM probe. Concentrations tested to be effective for detection include, 100, 200, 300, 400, 500, 600, 700, 800, 900 nM and any ratio between 1:1, 1:2, 1:3, 1:4 of primer concentrations. The activation of Taq polymerase is achieved at 95°C for 5 minutes is followed by forty-four cycles of denaturation at 95°C for 15 seconds and annealing and elongation at 65° C for 30 seconds.

Example 5: Detection of MLV in formalin fixed paraffin embedded tissue samples.

[0100] Tumors from mice were removed and divided into 2 equal parts. One part of the tumor was formalin-fixed, paraffin-embedded and the other part of the tumor was frozen at -80°C. The FFPE mouse tumor tissue was cut in half, with one half spiked-in with a known copy number amount of pAZ3-emd and the other half was not spiked-in. A known copy number amount of pAZ3-emd was spiked-in to the frozen fresh mouse tissue and pre-processing incubation buffer. The FFPE and frozen fresh mouse tissues were incubated at 56°C overnight in a pre-processing incubation buffer containing proteinase K and dithiothreitol (DTT). The following day, the mouse tissue was processed on the Maxwell 16 instrument to extract out gDNA as per standard procedure. The extracted gDNA concentration was quantified on the

Nanodrop 1000. The extracted DNA was tested for presence of MLV and env2 sequences by qPCR with the results shown in Figure 14.

Example 6: XMRV/MLV RT-PCR Assay.

[0101] The XMRV (xenotropic murine leukemia virus-related virus) RT-PCR assay is performed to quantify RNA from biological samples, such as whole blood and plasma, without the need for RNA extraction. The assay employs a two-step amplification process with the initial step consisting of the distribution of 2 μ L of experimental sample directly into a cDNA reaction mix. Following completion of the reverse transcriptase (RT) cDNA synthesis, a 2 μ L aliquot is removed, transferred into a qPCR reaction mix and a qPCR protocol is performed. A quantitation curve is established with 7 non-zero samples comprising of serial dilutions of defined copy number of reference vector to generate a Ct value versus copy number correlation. Linear regression analysis generates an equation which is used to calculate the copy number in the sample. (Figure 15).

[0102] Four control samples and one reagent control are used for this assay and are run in parallel with all test samples. The two step reaction requires controls for both the RT procedure and the qPCR procedure. Therefore, a positive, a negative and a non-template control are included for the cDNA synthesis step and a positive, negative and non-template control are included for the qPCR portion of the process.

[0103] A negative matrix sample (*i.e.* whole blood) is spiked with a defined quantity of 22Rv1 viral vector (see description under 'Reference Standard'). This control is prepared fresh with each run to determine the efficiency of the cDNA generation in the RT step.

[0104] 22rev1 genomic DNA containing the integrated retroviral vector sequences of XMRV provides the best biological mimic of the actual amplification target to be screened in patient tissues.

[0105] A negative matrix sample (*i.e.* whole blood) is used as a negative RT control as it does not contain any detectable XMRV endogenous sequences. This control is prepared fresh with each run to verify that non-specific products are not generated during cDNA synthesis of the qRT step. Confirmation is obtained upon completion of the qPCR procedure. No amplification is expected.

[0106] DNA isolated from non-infected U-87 cell is used as a negative control as it does not contain any endogenous sequences detectable by the XMRV primer sets.

[0107] The 22Rv1 human prostate carcinoma epithelial cell line has been shown to produce high-titer of the human retrovirus XMRV. This cell line was bought from ATCC and propagated in RPMI-1640 Medium containing 10% FBS, Sodium Pyruvate and Glutamax. The cell line was passaged four times before obtaining the supernatant containing the viral vector. The supernatant was filtered through a 0.45 μ m filter and stored at -80°C.

[0108] Reference vector 22Rv1 was used to spike PBS for generating a quantitation curve. Known copy numbers of vector were serially diluted to generate a Ct value versus copy number correlation. Linear regression analysis generates an equation which was used to calculate the copy number in the sample. Copy number was determined by a titer analysis which measures the number of copies of the viral genome integrated into the genome of target cells (transduction units, TU). The copy number was measured in TU equivalents.

[0109] Several studies were conducted to determine the appropriate primer sets, the optimal concentration for the reactions and the optimal temperature for the cycling parameters. Specific primer sets were designed and tested for human derived material. The goal of these experiments was to identify primer sets that were XMRV specific and did not present background in test samples.

[0110] The following primer sets were identified for targeting genes specific for XMRV:

1. XMRV gag

XMRV 628F (5'-ACTACCCCTCTGAGTCTAACC-3') (SEQ ID NO:1)
 XMRV764R (5'-GGCCATCCTACATTGAAAGTTG-3') (SEQ ID NO:2)
 XMRV gag probe (5'-**FAM**-CGCATTGCATCCAACCACTCTGTG-3'-BHQ) (SEQ ID NO:3)

2. XMRV env

XMRV 6252F (5'-TTTGATTCTCAGTGGGCTC-3') (SEQ ID NO:4)
 XMRV6391R (5'-CGATACAGTCTTAGTCCCATG-3') (SEQ ID NO:5)
 XMRV env probe (5'-**HEX**-CCCTTTTACCCGCGTCAGTGAATTCT-3'-BHQ) (SEQ ID NO:6)

[0111] Two primer probe sets were designed for detection of all MLV related retroviruses and XMRV.

3. MLV Pol 1

EP 2 913 402 B1

pol-F (5'-AACAAGCGGGTGAAGACATC-3') (SEQ ID NO:7)
pol-R (5'-CAAAGGCGAAGAGAGGCTGAC-3') (SEQ ID NO:8)
pol probe (5'-HEX-CCCACCGTGCCCAACCTTACAACC-3'-TAMRA) (SEQ ID NO:9)

4. MLV Pol 2

5' Pol2 Primer(CAAGGGGCTACTGGAGGAAAG) (SEQ ID NO:10)
3' Pol2 Primer:(CTTTCCTCCATGTACCAGACTG) (SEQ ID NO:11)
Pol2 Probe:(5HEX/TATCGCTGGACCACGGATCGCAA/3BHQ_1) (SEQ ID NO:12)

[0112] Two primer probe sets were designed for detection of amphotropic MLV virus.

5. MLV Env2

5'Env2 primer: 5'-ACCCTCAACCGCCCCCTACAAGT-3'(SEQ ID NO: 13)
3'Env2 primer: 5'-GTTAAGCGCCTGATAGGCTC-3' (SEQ ID NO: 14)
Env2 probe: 5'-/FAM/CCCCAAATGAAAGACCCCCGCTGACG/BHQ/-3'(SEQ ID NO: 15)

6. MLV LTR

One primer probe set was designed for detection of MLV in the LTR sequence

5'MLVLTR primer: AGC CCA CAA CCC CTC ACT C (SEQ ID NO:16)
3' MLVLTR primer TCT CCC GAT CCC GGA CGA (SEQ ID NO:17)
MLVLTR probe: FAM-CCCCAAATGAAAGACCCCCGCTGACG 3BHQ_1 (SEQ ID NO:18)

7. Cytosine deaminase gene

One primer probe set was designed for detection of the Cytosine deaminase gene

5' yCD2 Primer: (ATC ATC ATG TAC GGC ATC CCT AG) (SEQ ID NO:19)
3' yCD2 Primer: (TGA ACT GCT TCA TCA GCT TCT TAC) (SEQ ID NO:20)
yCD2 Probe: (5FAM/TCA TCG TCA ACA ACC ACC ACC TCG T/3BHQ_1) (SEQ ID NO:21).

Example 7: Monitoring of GBM patients treated with an MLV vector.

[0113] Open-label, ascending-dose trial of the safety and tolerability of increasing doses of Toca 511 administered to subjects with recurrent High Grade Glioma (including GBM) who have undergone surgery followed by adjuvant radiation and chemotherapy was carried out (see <http://clinicaltrials.gov/ct2/show/NCT01156584?term=tocagen&rank=1>). Ascending doses of Toca 511 (aka T5.0002) were prepared suitable for clinical use (WO2010148203) and delivered via stereotactic transcranial injection into the tumor. The starting dose was 2.6×10^3 TU/g. Subjects meeting all of the inclusion and none of the exclusion criteria received Toca 511 via stereotactic, transcranial injection into their tumor. Approximately three weeks (± 1 week) later subjects underwent a baseline gadolinium-enhanced MRI (Gd-MRI) scan and then began treatment with oral 5-FC at approximately 130 mg/kg/day for 6 days. On the 4th, 5th or 6th day of dosing the trough 5-FC serum concentration was determined and the dose of 5-FC adjusted in subsequent cycles to maintain the trough concentration in the therapeutic range. If tolerated, these 6-day courses of 5-FC were repeated approximately every 4 weeks (± 1 week) until institution of new antineoplastic treatment for tumor progression. Subjects undergo Gd-MRI scanning approximately every 8 weeks. Tumor response are assessed using the Macdonald criteria. A standard dose-escalation algorithm is being followed. Three subjects are evaluated at each of up to four dose levels of Toca 511 (2.6×10^3 , 9.5×10^3 , 2.5×10^4 , and the Maximum Feasible Dose [MFD], not to exceed 1×10^5 TU/g). So far three patients at the lowest dose level have been treated. Two patients 101 and 102 were monitored using qPCR testing of whole blood DNA (MLVLTR Primer probe set) and RT-qPCR using the MLV env2 primer-probe set. In addition, saliva and urine were monitored by DNA qPCR, and antibodies to the vector were measured (ref MLV ELISA Application?). These data are shown in Figure 16.

[0114] Other primers useful in the methods and composition of the disclosure for detecting XMRV and MLV related viruses include those in Table 1.

Table 1:

Sequence Definition	Pair Rating	Product Length	Product T _m	Sense Primer	Anti-sense Primer	5 XMRV position	3 XMRV position
XMRV 1581-1778	66.8	178	78.2	AGGTAGGAACCACTAGTCC	AGGGTCATAAGGAGTGATCC	1581	1758
XMRV 1581-1778	63.6	168	78.3	AGGTAGGAACCACTAGTCC	GGAGTGATACCTGCGATAGGC	1581	1748
XMRV 1581-1778	63.5	198	78.9	AGGTAGGAACCACTAGTCC	GTTTCTTGCCCTGGGTCCTC	1581	1778
XMRV 1581-1778	62.5	188	78.8	AGGTAGGAACCACTAGTCC	CTGGGTCTCAGGGTCATAA	1581	1768
XMRV 1581-1778	60.1	173	78	AGGTAGGAACCACTAGTCC	CATAAGGAGTGATACCTGCGA	1581	1753
XMRV 1729-1948	73.3	195	75.8	TCGCAGGTACACTCCTTATG	TCTCTTTCTTCCGGGTTTC	1734	1928
XMRV 1729-1948	69.8	215	76	TCGCAGGTACACTCCTTATG	TTTCTCTCCTGATACGTTCC	1734	1948
XMRV 1729-1948	69	200	76	TCGCAGGTACACTCCTTATG	GTTCCCTCTCTTTCTTCCGGG	1734	1933
XMRV 1729-1948	68.2	200	76	GCCTATCGCAGGTACACTCC	TCTCTTTCTTCCGGGTTTC	1729	1928
XMRV 1729-1948	66	205	76.2	GCCTATCGCAGGTACACTCC	GTTCCCTCTCTTTCTTCCGGG	1729	1933
XMRV 1729-1948	66	210	76.2	TCGCAGGTACACTCCTTATG	CTCCTGATACGTTCTCTCT	1734	1943
XMRV 1729-1948	63.8	180	75.1	TTATGACCCTGAGGACCCAG	TCTCTTTCTTCCGGGTTTC	1749	1928
XMRV 1729-1948	61.8	200	75.4	TTATGACCCTGAGGACCCAG	TTTCTCTCCTGATACGTTCC	1749	1948
XMRV 1729-1948	61.4	210	75.8	GGTACACTCCTTATGACCCT	TTTCTCTCCTGATACGTTCC	1739	1948
XMRV 1729-1948	61.1	185	75.3	TTATGACCCTGAGGACCCAG	GTTCCCTCTCTTTCTTCCGGG	1749	1933

EP 2 913 402 B1

[0115] Table 2 provides additional primer pairs and probes that are useful in the methods and compositions of the disclosure, for example for quantitative PCR and quantitative RT-PCR, and useful for XMRV or XMRV *and* MLV related virus detection:

5

10

15

20

25

30

35

40

45

50

55

Table 2:

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5'XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	82.7	TAGCGGCGGACCTCTCATTGACCTT	66.7	56	106	CTATAAAGTCCAAACCTCCTAAGCC	GAGGAAGGTTGTGCTCCGTAC	1015	1121	no
XMRV gag	84.9	AGCGGCGGACCTCTCATTGACCTTC	67.6	60	114	TCTATAAAGTCCAAACCTCCTAAGC	GGCAGAGGAGGAGAGTTGTG	1015	1128	no
XMRV gag	83.5	AGCGGCGGACCTCTCATTGACCTT	67	58.3	114	TCTATAAAGTCCAAACCTCCTAAGC	GGCAGAGGAGGAGAGTTGTG	1015	1128	no
XMRV gag	82.7	TAGCGGCGGACCTCTCATTGACCTT	66.7	56	114	TCTATAAAGTCCAAACCTCCTAAGC	GGCAGAGGAGGAGAGTTGTG	1015	1128	no
XMRV gag	82.5	AGCGGCGGACCTCTCATTGACCT	66.6	60.9	114	TCTATAAAGTCCAAACCTCCTAAGC	GGCAGAGGAGGAGAGTTGTG	1015	1128	no
XMRV gag	82.3	ATAGCGGCGGACCTTCTCATTGACCT	66.5	56	114	TCTATAAAGTCCAAACCTCCTAAGC	GGCAGAGGAGGAGAGTTGTG	1015	1128	no
XMRV gag	81.6	TAGCGGCGGACCTCTCATTGACCT	66.3	58.3	114	TCTATAAAGTCCAAACCTCCTAAGC	GGCAGAGGAGGAGAGTTGTG	1015	1128	no
XMRV gag	84.9	AGCGGCGGACCTCTCATTGACCTTC	67.6	60	129	TCTATAAAGTCCAAACCTCCTAAGC	TTCATTGTTCTCTCCCTGGCAGAG	1015	1143	no
XMRV gag	83.5	AGCGGCGGACCTCTCATTGACCTT	67	58.3	129	TCTATAAAGTCCAAACCTCCTAAGC	TTCATTGTTCTCTCCCTGGCAGAG	1015	1143	no
XMRV gag	82.7	TAGCGGCGGACCTCTCATTGACCTT	66.7	56	129	TCTATAAAGTCCAAACCTCCTAAGC	TTCATTGTTCTCTCCCTGGCAGAG	1015	1143	no
XMRV gag	82.5	AGCGGCGGACCTCTCATTGACCT	66.6	60.9	129	TCTATAAAGTCCAAACCTCCTAAGC	TTCATTGTTCTCTCCCTGGCAGAG	1015	1143	no
XMRV gag	82.3	ATAGCGGCGGACCTTCTCATTGACCT	66.5	56	129	TCTATAAAGTCCAAACCTCCTAAGC	TTCATTGTTCTCTCCCTGGCAGAG	1015	1143	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	81.6	TAGCGGCGGACCTC CTCATTGACCT	66.3	58.3	129	TCTATAAAGTCCAAA CCTCCTAAGC	TTCATTGTTCTCCC TGCGAGAG	1015	1143	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	139	TCTATAAAGTCCAAA CCTCCTAAGC	CCGCCTCTTCTTCA TTGTTCTC	1015	1153	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	139	TCTATAAAGTCCAAA CCTCCTAAGC	CCGCCTCTTCTTCA TTGTTCTC	1015	1153	no
XMRV gag	82.7	TAGCGGCGGACCTC CTCATTGACCTT	66.7	56	139	TCTATAAAGTCCAAA CCTCCTAAGC	CCGCCTCTTCTTCA TTGTTCTC	1015	1153	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	139	TCTATAAAGTCCAAA CCTCCTAAGC	CCGCCTCTTCTTCA TTGTTCTC	1015	1153	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	139	TCTATAAAGTCCAAA CCTCCTAAGC	CCGCCTCTTCTTCA TTGTTCTC	1015	1153	no
XMRV gag	81.6	TAGCGGCGGACCTC CTCATTGACCT	66.3	58.3	139	TCTATAAAGTCCAAA CCTCCTAAGC	CCGCCTCTTCTTCA TTGTTCTC	1015	1153	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	140	TCTATAAAGTCCAAA CCTCCTAAGC	GCCGCCTCTTCTTC ATTGTTTC	1015	1154	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	140	TCTATAAAGTCCAAA CCTCCTAAGC	GCCGCCTCTTCTTC ATTGTTTC	1015	1154	no
XMRV gag	82.7	TAGCGGCGGACCTC CTCATTGACCTT	66.7	56	140	TCTATAAAGTCCAAA CCTCCTAAGC	GCCGCCTCTTCTTC ATTGTTTC	1015	1154	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	140	TCTATAAAGTCCAAA CCTCCTAAGC	GCCGCCTCTTCTTC ATTGTTTC	1015	1154	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	140	TCTATAAAGTCCAAA CCTCCTAAGC	GCCGCCTCTTCTTC ATTGTTTC	1015	1154	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	81.6	TAGCGGGGACCTCTCATTGACCT	66.3	58.3	140	TCTATAAAGTCCAAACCTCCTAAGC	GCCGCCTCTTCTTCATTG	1015	1154	no
XMRV gag	84.9	AGCGGGGACCTCTCATTGACCTTC	67.6	60	142	TCTATAAAGTCCAAACCTCCTAAGC	TGCGCGCCTCTTCTTCATTG	1015	1156	no
XMRV gag	83.5	AGCGGGGACCTCTCATTGACCTT	67	58.3	142	TCTATAAAGTCCAAACCTCCTAAGC	TGCGCGCCTCTTCTTCATTG	1015	1156	no
XMRV gag	82.7	TAGCGGGGACCTCTCATTGACCTT	66.7	56	142	TCTATAAAGTCCAAACCTCCTAAGC	TGCGCGCCTCTTCTTCATTG	1015	1156	no
XMRV gag	82.5	AGCGGGGACCTCTCATTGACCT	66.6	60.9	142	TCTATAAAGTCCAAACCTCCTAAGC	TGCGCGCCTCTTCTTCATTG	1015	1156	no
XMRV gag	82.3	ATAGCGGGGACCTCTCATTGACCT	66.5	56	142	TCTATAAAGTCCAAACCTCCTAAGC	TGCGCGCCTCTTCTTCATTG	1015	1156	no
XMRV gag	81.6	TAGCGGGGACCTCTCATTGACCT	66.3	58.3	142	TCTATAAAGTCCAAACCTCCTAAGC	TGCGCGCCTCTTCTTCATTG	1015	1156	no
XMRV gag	84.9	AGCGGGGACCTCTCATTGACCTTC	67.6	60	144	TCTATAAAGTCCAAACCTCCTAAGC	GGTGGCCGCCTCTTCTTC	1015	1158	no
XMRV gag	83.5	AGCGGGGACCTCTCATTGACCTT	67	58.3	144	TCTATAAAGTCCAAACCTCCTAAGC	GGTGGCCGCCTCTTCTTC	1015	1158	no
XMRV gag	82.7	TAGCGGGGACCTCTCATTGACCTT	66.7	56	144	TCTATAAAGTCCAAACCTCCTAAGC	GGTGGCCGCCTCTTCTTC	1015	1158	no
XMRV gag	82.5	AGCGGGGACCTCTCATTGACCT	66.6	60.9	144	TCTATAAAGTCCAAACCTCCTAAGC	GGTGGCCGCCTCTTCTTC	1015	1158	no
XMRV gag	82.3	ATAGCGGGGACCTCTCATTGACCT	66.5	56	144	TCTATAAAGTCCAAACCTCCTAAGC	GGTGGCCGCCTCTTCTTC	1015	1158	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	81.6	TAGCGGGGACCTCTCATTGACCT	66.3	58.3	144	TCTATAAAGTCCAAACCTCCTAAGC	GGTGGCCGCCTCTTCTC	1015	1158	no
XMRV gag	82.7	TAGCGGGGACCTCTCATTGACCTT	66.7	56	188	CTATAAAGTCCAAACCTCCTAAGCC	CCGAGTCGAGACACCATG	1015	1203	no
XMRV gag	84.9	AGCGGGGACCTCTCATTGACCTTC	67.6	60	191	CTATAAAGTCCAAACCTCCTAAGCC	TCCCCGCAGTCGAGAC	1015	1206	no
XMRV gag	82.7	TAGCGGGGACCTCTCATTGACCTT	66.7	56	191	CTATAAAGTCCAAACCTCCTAAGCC	TCCCCGCAGTCGAGAC	1015	1206	no
XMRV gag	84.9	AGCGGGGACCTCTCATTGACCTTC	67.6	60	194	TCTATAAAGTCCAAACCTCCTAAGC	CTTCCCCGAGTCGAGAC	1015	1208	no
XMRV gag	83.5	AGCGGGGACCTCTCATTGACCTT	67	58.3	194	TCTATAAAGTCCAAACCTCCTAAGC	CTTCCCCGAGTCGAGAC	1015	1208	no
XMRV gag	82.7	TAGCGGGGACCTCTCATTGACCTT	66.7	56	194	TCTATAAAGTCCAAACCTCCTAAGC	CTTCCCCGAGTCGAGAC	1015	1208	no
XMRV gag	82.5	AGCGGGGACCTCTCATTGACCT	66.6	60.9	194	TCTATAAAGTCCAAACCTCCTAAGC	CTTCCCCGAGTCGAGAC	1015	1208	no
XMRV gag	82.3	ATAGCGGGGACCTCTCATTGACCT	66.5	56	194	TCTATAAAGTCCAAACCTCCTAAGC	CTTCCCCGAGTCGAGAC	1015	1208	no
XMRV gag	81.6	TAGCGGGGACCTCTCATTGACCT	66.3	58.3	194	TCTATAAAGTCCAAACCTCCTAAGC	CTTCCCCGAGTCGAGAC	1015	1208	no
XMRV gag	84.9	AGCGGGGACCTCTCATTGACCTTC	67.6	60	199	CTATAAAGTCCAAACCTCCTAAGCC	TCTCTCTTCCCCCGAGTC	1015	1214	no
XMRV gag	82.7	TAGCGGGGACCTCTCATTGACCTT	66.7	56	199	CTATAAAGTCCAAACCTCCTAAGCC	TCTCTCTTCCCCCGAGTC	1015	1214	no
XMRV gag	81.6	TAGCGGGGACCTCTCATTGACCT	66.3	58.3	199	CTATAAAGTCCAAACCTCCTAAGCC	TCTCTCTTCCCCCGAGTC	1015	1214	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
		CTCATTGACCT				CTCCTAAGCC	CAGTC			
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	106	CTATAAAGTCCAAAC CTCCTAAGCC	GAGGAAGGTTGT GCTCCGTAC	1016	1121	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	106	CTATAAAGTCCAAAC CTCCTAAGCC	GAGGAAGGTTGT GCTCCGTAC	1016	1121	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	106	CTATAAAGTCCAAAC CTCCTAAGCC	GAGGAAGGTTGT GCTCCGTAC	1016	1121	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	106	CTATAAAGTCCAAAC CTCCTAAGCC	GAGGAAGGTTGT GCTCCGTAC	1016	1121	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	106	CTATAAAGTCCAAAC CTCCTAAGCC	GAGGAAGGTTGT GCTCCGTAC	1016	1121	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	111	CTATAAAGTCCAAAC CTCCTAAGCC	CAGAGGAGGAAG GTTGTGCTC	1016	1126	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	111	CTATAAAGTCCAAAC CTCCTAAGCC	CAGAGGAGGAAG GTTGTGCTC	1016	1126	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	111	CTATAAAGTCCAAAC CTCCTAAGCC	CAGAGGAGGAAG GTTGTGCTC	1016	1126	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	111	CTATAAAGTCCAAAC CTCCTAAGCC	CAGAGGAGGAAG GTTGTGCTC	1016	1126	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	111	CTATAAAGTCCAAAC CTCCTAAGCC	CAGAGGAGGAAG GTTGTGCTC	1016	1126	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	111	CTATAAAGTCCAAAC CTCCTAAGCC	CAGAGGAGGAAG GTTGTGCTC	1016	1126	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	188	CTATAAAGTCCAAAC CTCCTAAGCC	CCGCAGTCGAGAC ACCATG	1016	1203	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	188	CTATAAAGTCCAAAC CTCCTAAGCC	CCGCAGTCGAGAC ACCATG	1016	1203	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	188	CTATAAAGTCCAAAC CTCCTAAGCC	CCGCAGTCGAGAC ACCATG	1016	1203	no
XMRV gag	82.3	ATAGCGGCGGACCTC TCTCATTGACCT	66.5	56	188	CTATAAAGTCCAAAC CTCCTAAGCC	CCGCAGTCGAGAC ACCATG	1016	1203	no
XMRV gag	81.6	TAGCGGCGGACCTC CTCATTGACCT	66.3	58.3	188	CTATAAAGTCCAAAC CTCCTAAGCC	CCGCAGTCGAGAC ACCATG	1016	1203	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	191	CTATAAAGTCCAAAC CTCCTAAGCC	TCCCCGCAGTCGA GACAC	1016	1206	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	191	CTATAAAGTCCAAAC CTCCTAAGCC	TCCCCGCAGTCGA GACAC	1016	1206	no
XMRV gag	82.3	ATAGCGGCGGACCTC TCTCATTGACCT	66.5	56	191	CTATAAAGTCCAAAC CTCCTAAGCC	TCCCCGCAGTCGA GACAC	1016	1206	no
XMRV gag	81.6	TAGCGGCGGACCTC CTCATTGACCT	66.3	58.3	191	CTATAAAGTCCAAAC CTCCTAAGCC	TCCCCGCAGTCGA GACAC	1016	1206	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTT	67.6	60	195	CTATAAAGTCCAAAC CTCCTAAGCC	TCCTTCCCCGCAGT CGAG	1016	1210	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	195	CTATAAAGTCCAAAC CTCCTAAGCC	TCCTTCCCCGCAGT CGAG	1016	1210	no
XMRV gag	82.7	TAGCGGCGGACCTC CTCATTGACCTT	66.7	56	195	CTATAAAGTCCAAAC CTCCTAAGCC	TCCTTCCCCGCAGT CGAG	1016	1210	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	195	CTATAAAGTCCAAAC CTCCTAAGCC	TCCTTCCCCGCAGT CGAG	1016	1210	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	82.3	ATAGCGGGGACCTCTCATTTGACCT	66.5	56	195	CTATAAAGTCCAAACCTCCTAAGCC	TCCTTCCCCGCGAGT	1016	1210	no
XMRV gag	81.6	TAGCGGGGACCTCTCATTTGACCT	66.3	58.3	195	CTATAAAGTCCAAACCTCCTAAGCC	TCCTTCCCCGCGAGT	1016	1210	no
XMRV gag	83.5	AGCGGGGACCTCTCATTTGACCT	67	58.3	199	CTATAAAGTCCAAACCTCCTAAGCC	TCTCTCCTTCCCCGAGTC	1016	1214	no
XMRV gag	82.5	AGCGGGGACCTCTCATTTGACCT	66.6	60.9	199	CTATAAAGTCCAAACCTCCTAAGCC	TCTCTCCTTCCCCGAGTC	1016	1214	no
XMRV gag	82.3	ATAGCGGGGACCTCTCATTTGACCT	66.5	56	199	CTATAAAGTCCAAACCTCCTAAGCC	TCTCTCCTTCCCCGAGTC	1016	1214	no
XMRV gag	84.9	AGCGGGGACCTCTCATTTGACCT	67.6	60	127	TATAAAGTCCAAACCTCCTAAGCC	TTCATTGTTCTCCC	1017	1143	no
XMRV gag	83.5	AGCGGGGACCTCTCATTTGACCT	67	58.3	127	TATAAAGTCCAAACCTCCTAAGCC	TTCATTGTTCTCCC	1017	1143	no
XMRV gag	82.7	TAGCGGGGACCTCTCATTTGACCT	66.7	56	127	TATAAAGTCCAAACCTCCTAAGCC	TTCATTGTTCTCCC	1017	1143	no
XMRV gag	82.5	AGCGGGGACCTCTCATTTGACCT	66.6	60.9	127	TATAAAGTCCAAACCTCCTAAGCC	TTCATTGTTCTCCC	1017	1143	no
XMRV gag	82.3	ATAGCGGGGACCTCTCATTTGACCT	66.5	56	127	TATAAAGTCCAAACCTCCTAAGCC	TTCATTGTTCTCCC	1017	1143	no
XMRV gag	81.6	TAGCGGGGACCTCTCATTTGACCT	66.3	58.3	127	TATAAAGTCCAAACCTCCTAAGCC	TTCATTGTTCTCCC	1017	1143	no
XMRV gag	84.9	AGCGGGGACCTCTCATTTGACCT	67.6	60	137	TATAAAGTCCAAACCTCCTAAGCC	CCGCCTCTTCTTCA	1017	1153	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	137	TATAAAGTCCAAACC TCCTAAGCC	CCGCCTCTTCTTCA TTGTTCTC	1017	1153	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	137	TATAAAGTCCAAACC TCCTAAGCC	CCGCCTCTTCTTCA TTGTTCTC	1017	1153	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	137	TATAAAGTCCAAACC TCCTAAGCC	CCGCCTCTTCTTCA TTGTTCTC	1017	1153	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	137	TATAAAGTCCAAACC TCCTAAGCC	CCGCCTCTTCTTCA TTGTTCTC	1017	1153	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	137	TATAAAGTCCAAACC TCCTAAGCC	CCGCCTCTTCTTCA TTGTTCTC	1017	1153	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	138	TATAAAGTCCAAACC TCCTAAGCC	GCCGCCTCTTCTTC ATTGTTTC	1017	1154	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	138	TATAAAGTCCAAACC TCCTAAGCC	GCCGCCTCTTCTTC ATTGTTTC	1017	1154	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	138	TATAAAGTCCAAACC TCCTAAGCC	GCCGCCTCTTCTTC ATTGTTTC	1017	1154	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	138	TATAAAGTCCAAACC TCCTAAGCC	GCCGCCTCTTCTTC ATTGTTTC	1017	1154	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	138	TATAAAGTCCAAACC TCCTAAGCC	GCCGCCTCTTCTTC ATTGTTTC	1017	1154	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	138	TATAAAGTCCAAACC TCCTAAGCC	GCCGCCTCTTCTTC ATTGTTTC	1017	1154	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	140	TATAAAGTCCAAACC TCCTAAGCC	TGCGCGCCTCTTCT TCATTG	1017	1156	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	140	TATAAAGTCCAAACC TCCTAAGCC	TGCGCGCCTCTTCT TCATTG	1017	1156	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	140	TATAAAGTCCAAACC TCCTAAGCC	TGCGCGCCTCTTCT TCATTG	1017	1156	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	140	TATAAAGTCCAAACC TCCTAAGCC	TGCGCGCCTCTTCT TCATTG	1017	1156	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	140	TATAAAGTCCAAACC TCCTAAGCC	TGCGCGCCTCTTCT TCATTG	1017	1156	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	140	TATAAAGTCCAAACC TCCTAAGCC	TGCGCGCCTCTTCT TCATTG	1017	1156	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	142	TATAAAGTCCAAACC TCCTAAGCC	GGTGGCCGCCTCT TCTTC	1017	1158	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	142	TATAAAGTCCAAACC TCCTAAGCC	GGTGGCCGCCTCT TCTTC	1017	1158	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	142	TATAAAGTCCAAACC TCCTAAGCC	GGTGGCCGCCTCT TCTTC	1017	1158	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	142	TATAAAGTCCAAACC TCCTAAGCC	GGTGGCCGCCTCT TCTTC	1017	1158	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	142	TATAAAGTCCAAACC TCCTAAGCC	GGTGGCCGCCTCT TCTTC	1017	1158	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	142	TATAAAGTCCAAACC TCCTAAGCC	GGTGGCCGCCTCT TCTTC	1017	1158	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	186	ATAAAGTCCAAACCT CCTAAGCC	CCGCAGTCGAGAC ACCATG	1018	1203	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	186	ATAAAGTCCAAACCT CCTAAGCC	CCGAGTCGAGAC ACCATG	1018	1203	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	186	ATAAAGTCCAAACCT CCTAAGCC	CCGAGTCGAGAC ACCATG	1018	1203	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	186	ATAAAGTCCAAACCT CCTAAGCC	CCGAGTCGAGAC ACCATG	1018	1203	no
XMRV gag	82.3	ATAGCGGCGGAC TCTCATTGACCT	66.5	56	186	ATAAAGTCCAAACCT CCTAAGCC	CCGAGTCGAGAC ACCATG	1018	1203	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	186	ATAAAGTCCAAACCT CCTAAGCC	CCGAGTCGAGAC ACCATG	1018	1203	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	197	ATAAAGTCCAAACCT CCTAAGCC	TCTCTCCTTCCCCG CAGTC	1018	1214	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	197	ATAAAGTCCAAACCT CCTAAGCC	TCTCTCCTTCCCCG CAGTC	1018	1214	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	197	ATAAAGTCCAAACCT CCTAAGCC	TCTCTCCTTCCCCG CAGTC	1018	1214	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	197	ATAAAGTCCAAACCT CCTAAGCC	TCTCTCCTTCCCCG CAGTC	1018	1214	no
XMRV gag	82.3	ATAGCGGCGGAC TCTCATTGACCT	66.5	56	197	ATAAAGTCCAAACCT CCTAAGCC	TCTCTCCTTCCCCG CAGTC	1018	1214	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	197	ATAAAGTCCAAACCT CCTAAGCC	TCTCTCCTTCCCCG CAGTC	1018	1214	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	135	TAAAGTCCAAACCTC CTAAGCC	CCGCCTCTTCTTCA TTGTCTC	1019	1153	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGGGACCTC TCATTGACCTT	67	58.3	135	TAAAGTCCAAACCTC CTAAGGCC	CCGCCTCTTCTTCA TTGTTCTC	1019	1153	no
XMRV gag	82.7	TAGCGGGGACCT CTCATTGACCTT	66.7	56	135	TAAAGTCCAAACCTC CTAAGGCC	CCGCCTCTTCTTCA TTGTTCTC	1019	1153	no
XMRV gag	82.5	AGCGGGGACCTC TCATTGACCT	66.6	60.9	135	TAAAGTCCAAACCTC CTAAGGCC	CCGCCTCTTCTTCA TTGTTCTC	1019	1153	no
XMRV gag	82.3	ATAGCGGGGACCT TCTCATTGACCT	66.5	56	135	TAAAGTCCAAACCTC CTAAGGCC	CCGCCTCTTCTTCA TTGTTCTC	1019	1153	no
XMRV gag	81.6	TAGCGGGGACCT CTCATTGACCT	66.3	58.3	135	TAAAGTCCAAACCTC CTAAGGCC	CCGCCTCTTCTTCA TTGTTCTC	1019	1153	no
XMRV gag	84.9	AGCGGGGACCTC TCATTGACCTTC	67.6	60	136	TAAAGTCCAAACCTC CTAAGGCC	GCCGCCTCTTCTTC ATTGTTTC	1019	1154	no
XMRV gag	83.5	AGCGGGGACCTC TCATTGACCTT	67	58.3	136	TAAAGTCCAAACCTC CTAAGGCC	GCCGCCTCTTCTTC ATTGTTTC	1019	1154	no
XMRV gag	82.7	TAGCGGGGACCT CTCATTGACCTT	66.7	56	136	TAAAGTCCAAACCTC CTAAGGCC	GCCGCCTCTTCTTC ATTGTTTC	1019	1154	no
XMRV gag	82.5	AGCGGGGACCTC TCATTGACCT	66.6	60.9	136	TAAAGTCCAAACCTC CTAAGGCC	GCCGCCTCTTCTTC ATTGTTTC	1019	1154	no
XMRV gag	81.6	TAGCGGGGACCT CTCATTGACCT	66.3	58.3	136	TAAAGTCCAAACCTC CTAAGGCC	GCCGCCTCTTCTTC ATTGTTTC	1019	1154	no
XMRV gag	84.9	AGCGGGGACCTC TCATTGACCTTC	67.6	60	140	TAAAGTCCAAACCTC CTAAGGCC	GGTGGCCGCCTCT TCTTC	1019	1158	no
XMRV gag	83.5	AGCGGGGACCTC TCATTGACCTT	67	58.3	140	TAAAGTCCAAACCTC CTAAGGCC	GGTGGCCGCCTCT TCTTC	1019	1158	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	82.7	TAGCGGCGGACCTCTCATTGACCTT	66.7	56	140	TAAAGTCCAAACCTCCTAAGCC	GGTGGCCGCCTCTCTTC	1019	1158	no
XMRV gag	82.5	AGCGGCGGACCTCTCATTGACCT	66.6	60.9	140	TAAAGTCCAAACCTCCTAAGCC	GGTGGCCGCCTCTCTTC	1019	1158	no
XMRV gag	82.3	ATAGCGGCGGACCTTCTCATTGACCT	66.5	56	136	TAAAGTCCAAACCTCCTAAGCC	GCCGCCTCTTCTTTCATTGTTTC	1019	1158	no
XMRV gag	82.3	ATAGCGGCGGACCTTCTCATTGACCT	66.5	56	140	TAAAGTCCAAACCTCCTAAGCC	GGTGGCCGCCTCTCTTC	1019	1158	no
XMRV gag	81.6	TAGCGGCGGACCTCTCATTGACCT	66.3	58.3	140	TAAAGTCCAAACCTCCTAAGCC	GGTGGCCGCCTCTCTTC	1019	1158	no
XMRV gag	84.9	AGCGGCGGACCTCTCATTGACCTTC	67.6	60	92	CCTCCTAAGCCCCCAGGTTCTC	GAGGAAGGTTGTGCTCCGTAC	1030	1121	no
XMRV gag	83.5	AGCGGCGGACCTCTCATTGACCTT	67	58.3	92	CCTCCTAAGCCCCCAGGTTCTC	GAGGAAGGTTGTGCTCCGTAC	1030	1121	no
XMRV gag	82.7	TAGCGGCGGACCTCTCATTGACCTT	66.7	56	92	CCTCCTAAGCCCCCAGGTTCTC	GAGGAAGGTTGTGCTCCGTAC	1030	1121	no
XMRV gag	82.5	AGCGGCGGACCTCTCATTGACCT	66.6	60.9	92	CCTCCTAAGCCCCCAGGTTCTC	GAGGAAGGTTGTGCTCCGTAC	1030	1121	no
XMRV gag	82.3	ATAGCGGCGGACCTTCTCATTGACCT	66.5	56	92	CCTCCTAAGCCCCCAGGTTCTC	GAGGAAGGTTGTGCTCCGTAC	1030	1121	no
XMRV gag	81.6	TAGCGGCGGACCTCTCATTGACCT	66.3	58.3	92	CCTCCTAAGCCCCCAGGTTCTC	GAGGAAGGTTGTGCTCCGTAC	1030	1121	no
XMRV gag	84.9	AGCGGCGGACCTCTCATTGACCTTC	67.6	60	99	CCTCCTAAGCCCCCAGGTTCTC	GGCAGAGGAGGAGGAGGTTGTG	1030	1128	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	99	CCTCCTAAGCCCCCAG GTTCTC	GGCAGAGGAGGA AGGTTGTG	1030	1128	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	99	CCTCCTAAGCCCCCAG GTTCTC	GGCAGAGGAGGA AGGTTGTG	1030	1128	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	99	CCTCCTAAGCCCCCAG GTTCTC	GGCAGAGGAGGA AGGTTGTG	1030	1128	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	99	CCTCCTAAGCCCCCAG GTTCTC	GGCAGAGGAGGA AGGTTGTG	1030	1128	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	99	CCTCCTAAGCCCCCAG GTTCTC	GGCAGAGGAGGA AGGTTGTG	1030	1128	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	114	CCTCCTAAGCCCCCAG GTTCTC	TTCATTGTTCTCCC TGGCAGAG	1030	1143	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	114	CCTCCTAAGCCCCCAG GTTCTC	TTCATTGTTCTCCC TGGCAGAG	1030	1143	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	114	CCTCCTAAGCCCCCAG GTTCTC	TTCATTGTTCTCCC TGGCAGAG	1030	1143	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	114	CCTCCTAAGCCCCCAG GTTCTC	TTCATTGTTCTCCC TGGCAGAG	1030	1143	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	114	CCTCCTAAGCCCCCAG GTTCTC	TTCATTGTTCTCCC TGGCAGAG	1030	1143	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	124	CCTCCTAAGCCCCCAG GTTCTC	CCGCCTCTTCTTCA TTGTCTC	1030	1153	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	124	CCTCCTAAGCCCCCAG GTTCTC	CCGCCTCTTCTTCA TTGTCTC	1030	1153	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	82.7	TAGCGGCGGACCTCTCATTGACCTT	66.7	56	124	CCTCCTAAGCCCCCAGGTTCTC	CCGCCTCTTCTTCA TTGTTCTC	1030	1153	no
XMRV gag	82.5	AGCGGCGGACCTCTCATTGACCT	66.6	60.9	124	CCTCCTAAGCCCCCAGGTTCTC	CCGCCTCTTCTTCA TTGTTCTC	1030	1153	no
XMRV gag	82.3	ATAGCGGCGGACCTTCTCATTGACCT	66.5	56	124	CCTCCTAAGCCCCCAGGTTCTC	CCGCCTCTTCTTCA TTGTTCTC	1030	1153	no
XMRV gag	81.6	TAGCGGCGGACCTCTCATTGACCT	66.3	58.3	124	CCTCCTAAGCCCCCAGGTTCTC	CCGCCTCTTCTTCA TTGTTCTC	1030	1153	no
XMRV gag	84.9	AGCGGCGGACCTCTCATTGACCTTC	67.6	60	125	CCTCCTAAGCCCCCAGGTTCTC	GCCGCCTCTTCTTC ATTGTTTC	1030	1154	no
XMRV gag	83.5	AGCGGCGGACCTCTCATTGACCTT	67	58.3	125	CCTCCTAAGCCCCCAGGTTCTC	GCCGCCTCTTCTTC ATTGTTTC	1030	1154	no
XMRV gag	82.7	TAGCGGCGGACCTCTCATTGACCTT	66.7	56	125	CCTCCTAAGCCCCCAGGTTCTC	GCCGCCTCTTCTTC ATTGTTTC	1030	1154	no
XMRV gag	82.5	AGCGGCGGACCTCTCATTGACCT	66.6	60.9	125	CCTCCTAAGCCCCCAGGTTCTC	GCCGCCTCTTCTTC ATTGTTTC	1030	1154	no
XMRV gag	82.3	ATAGCGGCGGACCTTCTCATTGACCT	66.5	56	125	CCTCCTAAGCCCCCAGGTTCTC	GCCGCCTCTTCTTC ATTGTTTC	1030	1154	no
XMRV gag	81.6	TAGCGGCGGACCTCTCATTGACCT	66.3	58.3	125	CCTCCTAAGCCCCCAGGTTCTC	GCCGCCTCTTCTTC ATTGTTTC	1030	1154	no
XMRV gag	84.9	AGCGGCGGACCTCTCATTGACCTTC	67.6	60	127	CCTCCTAAGCCCCCAGGTTCTC	TGCCCGCCTCTTCT TCATTG	1030	1156	no
XMRV gag	83.5	AGCGGCGGACCTCTCATTGACCTT	67	58.3	127	CCTCCTAAGCCCCCAGGTTCTC	TGCCCGCCTCTTCT TCATTG	1030	1156	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	82.7	TAGCGGCGGACCTC TCATTGACCTT	66.7	56	127	CCTCCTAAGCCCCCAG GTTCTC	TGGCCGCCTCTTCT TCATTG	1030	1156	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	127	CCTCCTAAGCCCCCAG GTTCTC	TGGCCGCCTCTTCT TCATTG	1030	1156	no
XMRV gag	82.3	ATAGCGGCGGACCTC TCTCATTGACCT	66.5	56	127	CCTCCTAAGCCCCCAG GTTCTC	TGGCCGCCTCTTCT TCATTG	1030	1156	no
XMRV gag	81.6	TAGCGGCGGACCTC CTCATTGACCT	66.3	58.3	127	CCTCCTAAGCCCCCAG GTTCTC	TGGCCGCCTCTTCT TCATTG	1030	1156	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	132	CCTCCTAAGCCCCCAG GTTCTC	GGTGGTGCCCGCC TCTTC	1030	1158	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	129	CCTCCTAAGCCCCCAG GTTCTC	GGTGGCCGCCTCT TCTTC	1030	1158	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	132	CCTCCTAAGCCCCCAG GTTCTC	GGTGGTGCCCGCC TCTTC	1030	1158	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	129	CCTCCTAAGCCCCCAG GTTCTC	GGTGGCCGCCTCT TCTTC	1030	1158	no
XMRV gag	82.7	TAGCGGCGGACCTC CTCATTGACCTT	66.7	56	132	CCTCCTAAGCCCCCAG GTTCTC	GGTGGTGCCCGCC TCTTC	1030	1158	no
XMRV gag	82.7	TAGCGGCGGACCTC CTCATTGACCTT	66.7	56	129	CCTCCTAAGCCCCCAG GTTCTC	GGTGGCCGCCTCT TCTTC	1030	1158	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	132	CCTCCTAAGCCCCCAG GTTCTC	GGTGGTGCCCGCC TCTTC	1030	1158	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	129	CCTCCTAAGCCCCCAG GTTCTC	GGTGGCCGCCTCT TCTTC	1030	1158	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	82.3	ATAGCGGCGGACCTCTCATTGACCT	66.5	56	132	CCTCCTAAGCCCCCAGGTTCTC	GGTGGTGCCCGCC TCTTC	1030	1158	no
XMRV gag	82.3	ATAGCGGCGGACCTCTCATTGACCT	66.5	56	129	CCTCCTAAGCCCCCAGGTTCTC	GGTGGCGCCTCT TCTTC	1030	1158	no
XMRV gag	81.6	TAGCGGCGGACCTCTCATTGACCT	66.3	58.3	132	CCTCCTAAGCCCCCAGGTTCTC	GGTGGTGCCCGCC TCTTC	1030	1158	no
XMRV gag	81.6	TAGCGGCGGACCTCTCATTGACCT	66.3	58.3	129	CCTCCTAAGCCCCCAGGTTCTC	GGTGGCGCCTCT TCTTC	1030	1158	no
XMRV gag	84.9	AGCGGCGGACCTCTCATTGACCTTC	67.6	60	174	CCTCCTAAGCCCCCAGGTTCTC	CCGCAGTCGAGAC ACCATG	1030	1203	no
XMRV gag	83.5	AGCGGCGGACCTCTCATTGACCTT	67	58.3	174	CCTCCTAAGCCCCCAGGTTCTC	CCGCAGTCGAGAC ACCATG	1030	1203	no
XMRV gag	82.7	TAGCGGCGGACCTCTCATTGACCTT	66.7	56	174	CCTCCTAAGCCCCCAGGTTCTC	CCGCAGTCGAGAC ACCATG	1030	1203	no
XMRV gag	82.5	AGCGGCGGACCTCTCATTGACCT	66.6	60.9	174	CCTCCTAAGCCCCCAGGTTCTC	CCGCAGTCGAGAC ACCATG	1030	1203	no
XMRV gag	82.3	ATAGCGGCGGACCTCTCATTGACCT	66.5	56	174	CCTCCTAAGCCCCCAGGTTCTC	CCGCAGTCGAGAC ACCATG	1030	1203	no
XMRV gag	81.6	TAGCGGCGGACCTCTCATTGACCT	66.3	58.3	174	CCTCCTAAGCCCCCAGGTTCTC	CCGCAGTCGAGAC ACCATG	1030	1203	no
XMRV gag	84.9	AGCGGCGGACCTCTCATTGACCTTC	67.6	60	177	CCTCCTAAGCCCCCAGGTTCTC	TCCCCGACGTCGA GACAC	1030	1206	no
XMRV gag	83.5	AGCGGCGGACCTCTCATTGACCTT	67	58.3	177	CCTCCTAAGCCCCCAGGTTCTC	TCCCCGACGTCGA GACAC	1030	1206	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	82.7	TAGCGGGGACCTCTCATTGACCTT	66.7	56	177	CCTCCTAAGCCCCCAGGTTCTC	TCCCCGCAGTCGAGACAC	1030	1206	no
XMRV gag	82.5	AGCGGGGACCTCTCATTGACCT	66.6	60.9	177	CCTCCTAAGCCCCCAGGTTCTC	TCCCCGCAGTCGAGACAC	1030	1206	no
XMRV gag	82.3	ATAGCGGGGACCTCTCATTGACCT	66.5	56	177	CCTCCTAAGCCCCCAGGTTCTC	TCCCCGCAGTCGAGACAC	1030	1206	no
XMRV gag	81.6	TAGCGGGGACCTCTCATTGACCT	66.3	58.3	177	CCTCCTAAGCCCCCAGGTTCTC	TCCCCGCAGTCGAGACAC	1030	1206	no
XMRV gag	84.9	AGCGGGGACCTCTCATTGACCTTC	67.6	60	181	CCTCCTAAGCCCCCAGGTTCTC	TCCTTCCCCGCAGTCGAG	1030	1210	no
XMRV gag	83.5	AGCGGGGACCTCTCATTGACCTT	67	58.3	181	CCTCCTAAGCCCCCAGGTTCTC	TCCTTCCCCGCAGTCGAG	1030	1210	no
XMRV gag	82.7	TAGCGGGGACCTCTCATTGACCTT	66.7	56	181	CCTCCTAAGCCCCCAGGTTCTC	TCCTTCCCCGCAGTCGAG	1030	1210	no
XMRV gag	82.5	AGCGGGGACCTCTCATTGACCT	66.6	60.9	181	CCTCCTAAGCCCCCAGGTTCTC	TCCTTCCCCGCAGTCGAG	1030	1210	no
XMRV gag	82.3	ATAGCGGGGACCTCTCATTGACCT	66.5	56	181	CCTCCTAAGCCCCCAGGTTCTC	TCCTTCCCCGCAGTCGAG	1030	1210	no
XMRV gag	81.6	TAGCGGGGACCTCTCATTGACCT	66.3	58.3	114	CCTCCTAAGCCCCCAGGTTCTC	TTCATTGTTCTCCC TGCGACAG	1030	1210	no
XMRV gag	81.6	TAGCGGGGACCTCTCATTGACCT	66.3	58.3	181	CCTCCTAAGCCCCCAGGTTCTC	TCCTTCCCCGCAGTCGAG	1030	1210	no
XMRV gag	84.9	AGCGGGGACCTCTCATTGACCTTC	67.6	60	185	CCTCCTAAGCCCCCAGGTTCTC	TCTCTCTTCTCCCGCAGTC	1030	1214	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	185	CCTCCTAAGCCCCCAG GTTCTC	TCTCTCCTTCCCCG CAGTC	1030	1214	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	185	CCTCCTAAGCCCCCAG GTTCTC	TCTCTCCTTCCCCG CAGTC	1030	1214	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	185	CCTCCTAAGCCCCCAG GTTCTC	TCTCTCCTTCCCCG CAGTC	1030	1214	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	185	CCTCCTAAGCCCCCAG GTTCTC	TCTCTCCTTCCCCG CAGTC	1030	1214	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	185	CCTCCTAAGCCCCCAG GTTCTC	TCTCTCCTTCCCCG CAGTC	1030	1214	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	197	CCTCCTAAGCCCCCAG GTTCTC	GCTGCGGGAGGG TCTCTC	1030	1226	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	197	CCTCCTAAGCCCCCAG GTTCTC	GCTGCGGGAGGG TCTCTC	1030	1226	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	197	CCTCCTAAGCCCCCAG GTTCTC	GCTGCGGGAGGG TCTCTC	1030	1226	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	197	CCTCCTAAGCCCCCAG GTTCTC	GCTGCGGGAGGG TCTCTC	1030	1226	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	197	CCTCCTAAGCCCCCAG GTTCTC	GCTGCGGGAGGG TCTCTC	1030	1226	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	197	CCTCCTAAGCCCCCAG GTTCTC	GCTGCGGGAGGG TCTCTC	1030	1226	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	90	TCCTAAGCCCCCAGGT TCTCC	GAGGAAGGTTGT GCTCCGTAC	1032	1121	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	90	TCCTAAGCCCCCAGGT TCTCC	GAGGAAGGTTGT GCTCCGTAC	1032	1121	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	90	TCCTAAGCCCCCAGGT TCTCC	GAGGAAGGTTGT GCTCCGTAC	1032	1121	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	90	TCCTAAGCCCCCAGGT TCTCC	GAGGAAGGTTGT GCTCCGTAC	1032	1121	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	90	TCCTAAGCCCCCAGGT TCTCC	GAGGAAGGTTGT GCTCCGTAC	1032	1121	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	90	TCCTAAGCCCCCAGGT TCTCC	GAGGAAGGTTGT GCTCCGTAC	1032	1121	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	95	TCCTAAGCCCCCAGGT TCTCC	CAGAGGAGGAAG GTTGTGCTC	1032	1126	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	95	TCCTAAGCCCCCAGGT TCTCC	CAGAGGAGGAAG GTTGTGCTC	1032	1126	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	95	TCCTAAGCCCCCAGGT TCTCC	CAGAGGAGGAAG GTTGTGCTC	1032	1126	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	95	TCCTAAGCCCCCAGGT TCTCC	CAGAGGAGGAAG GTTGTGCTC	1032	1126	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	95	TCCTAAGCCCCCAGGT TCTCC	CAGAGGAGGAAG GTTGTGCTC	1032	1126	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	95	TCCTAAGCCCCCAGGT TCTCC	CAGAGGAGGAAG GTTGTGCTC	1032	1126	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	97	TCCTAAGCCCCCAGGT TCTCC	GGCAGAGGAGGA AGGTTGTG	1032	1128	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	97	TCCTAAGCCCCCAGGT TCTCC	GGCAGAGGAGGA AGGTTGTG	1032	1128	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	97	TCCTAAGCCCCCAGGT TCTCC	GGCAGAGGAGGA AGGTTGTG	1032	1128	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	97	TCCTAAGCCCCCAGGT TCTCC	GGCAGAGGAGGA AGGTTGTG	1032	1128	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	97	TCCTAAGCCCCCAGGT TCTCC	GGCAGAGGAGGA AGGTTGTG	1032	1128	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	97	TCCTAAGCCCCCAGGT TCTCC	GGCAGAGGAGGA AGGTTGTG	1032	1128	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	99	TCCTAAGCCCCCAGGT TCTCC	CTGGCAGAGGAG GAAGGTTG	1032	1130	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	99	TCCTAAGCCCCCAGGT TCTCC	CTGGCAGAGGAG GAAGGTTG	1032	1130	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	99	TCCTAAGCCCCCAGGT TCTCC	CTGGCAGAGGAG GAAGGTTG	1032	1130	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	99	TCCTAAGCCCCCAGGT TCTCC	CTGGCAGAGGAG GAAGGTTG	1032	1130	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	99	TCCTAAGCCCCCAGGT TCTCC	CTGGCAGAGGAG GAAGGTTG	1032	1130	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	99	TCCTAAGCCCCCAGGT TCTCC	CTGGCAGAGGAG GAAGGTTG	1032	1130	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	112	TCCTAAGCCCCCAGGT TCTCC	TTCATTGTTCTCCC TGGCAGAG	1032	1143	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	112	TCCTAAGCCCCCAGGT TCTCC	TTCATTGTTCTCCC TGGCAGAG	1032	1143	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	112	TCCTAAGCCCCCAGGT TCTCC	TTCATTGTTCTCCC TGGCAGAG	1032	1143	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	112	TCCTAAGCCCCCAGGT TCTCC	TTCATTGTTCTCCC TGGCAGAG	1032	1143	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	112	TCCTAAGCCCCCAGGT TCTCC	TTCATTGTTCTCCC TGGCAGAG	1032	1143	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	112	TCCTAAGCCCCCAGGT TCTCC	TTCATTGTTCTCCC TGGCAGAG	1032	1143	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	122	TCCTAAGCCCCCAGGT TCTCC	CCGCCTCTTCTTCA TTGTTCTC	1032	1153	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	122	TCCTAAGCCCCCAGGT TCTCC	CCGCCTCTTCTTCA TTGTTCTC	1032	1153	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	122	TCCTAAGCCCCCAGGT TCTCC	CCGCCTCTTCTTCA TTGTTCTC	1032	1153	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	122	TCCTAAGCCCCCAGGT TCTCC	CCGCCTCTTCTTCA TTGTTCTC	1032	1153	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	122	TCCTAAGCCCCCAGGT TCTCC	CCGCCTCTTCTTCA TTGTTCTC	1032	1153	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	122	TCCTAAGCCCCCAGGT TCTCC	CCGCCTCTTCTTCA TTGTTCTC	1032	1153	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	123	TCCTAAGCCCCCAGGT TCTCC	GCCGCCTCTTCTTC ATTGTTTC	1032	1154	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	123	TCCTAAGCCCCCAGGT TCTCC	GCCGCCTCTTCTTC ATTGTTT	1032	1154	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	123	TCCTAAGCCCCCAGGT TCTCC	GCCGCCTCTTCTTC ATTGTTT	1032	1154	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	123	TCCTAAGCCCCCAGGT TCTCC	GCCGCCTCTTCTTC ATTGTTT	1032	1154	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	123	TCCTAAGCCCCCAGGT TCTCC	GCCGCCTCTTCTTC ATTGTTT	1032	1154	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	123	TCCTAAGCCCCCAGGT TCTCC	GCCGCCTCTTCTTC ATTGTTT	1032	1154	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	125	TCCTAAGCCCCCAGGT TCTCC	TGGCCGCCTCTTCT TCATTG	1032	1156	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	172	TCCTAAGCCCCCAGGT TCTCC	CCGACGTCGAGAC ACCATG	1032	1156	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	125	TCCTAAGCCCCCAGGT TCTCC	TGGCCGCCTCTTCT TCATTG	1032	1156	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	125	TCCTAAGCCCCCAGGT TCTCC	TGGCCGCCTCTTCT TCATTG	1032	1156	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	125	TCCTAAGCCCCCAGGT TCTCC	TGGCCGCCTCTTCT TCATTG	1032	1156	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	125	TCCTAAGCCCCCAGGT TCTCC	TGGCCGCCTCTTCT TCATTG	1032	1156	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	125	TCCTAAGCCCCCAGGT TCTCC	TGGCCGCCTCTTCT TCATTG	1032	1156	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	127	TCCTAAGCCCCCAGGT TCTCC	GGTGGCCGCCTCT TCTTC	1032	1158	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	130	TCCTAAGCCCCCAGGT TCTCC	GGTGGTGGCCGCC TCTTC	1032	1158	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	127	TCCTAAGCCCCCAGGT TCTCC	GGTGGCCGCCTCT TCTTC	1032	1158	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	130	TCCTAAGCCCCCAGGT TCTCC	GGTGGTGGCCGCC TCTTC	1032	1158	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	127	TCCTAAGCCCCCAGGT TCTCC	GGTGGCCGCCTCT TCTTC	1032	1158	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	130	TCCTAAGCCCCCAGGT TCTCC	GGTGGTGGCCGCC TCTTC	1032	1158	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	127	TCCTAAGCCCCCAGGT TCTCC	GGTGGCCGCCTCT TCTTC	1032	1158	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	130	TCCTAAGCCCCCAGGT TCTCC	GGTGGTGGCCGCC TCTTC	1032	1158	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	127	TCCTAAGCCCCCAGGT TCTCC	GGTGGCCGCCTCT TCTTC	1032	1158	no
XMRV gag	82.3	ATAGCGGCGGACC TCTCATTGACCT	66.5	56	130	TCCTAAGCCCCCAGGT TCTCC	GGTGGTGGCCGCC TCTTC	1032	1158	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	127	TCCTAAGCCCCCAGGT TCTCC	GGTGGCCGCCTCT TCTTC	1032	1158	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	130	TCCTAAGCCCCCAGGT TCTCC	GGTGGTGGCCGCC TCTTC	1032	1161	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	172	TCCTAAGCCCCCAGGT TCTCC	CCGCAGTCGAGAC ACCATG	1032	1203	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	172	TCCTAAGCCCCCAGGT TCTCC	CCGCAGTCGAGAC ACCATG	1032	1203	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	172	TCCTAAGCCCCCAGGT TCTCC	CCGCAGTCGAGAC ACCATG	1032	1203	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	172	TCCTAAGCCCCCAGGT TCTCC	CCGCAGTCGAGAC ACCATG	1032	1203	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	172	TCCTAAGCCCCCAGGT TCTCC	CCGCAGTCGAGAC ACCATG	1032	1203	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	175	TCCTAAGCCCCCAGGT TCTCC	TCCCCGCGAGTCGA GACAC	1032	1206	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	175	TCCTAAGCCCCCAGGT TCTCC	TCCCCGCGAGTCGA GACAC	1032	1206	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	175	TCCTAAGCCCCCAGGT TCTCC	TCCCCGCGAGTCGA GACAC	1032	1206	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	175	TCCTAAGCCCCCAGGT TCTCC	TCCCCGCGAGTCGA GACAC	1032	1206	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	175	TCCTAAGCCCCCAGGT TCTCC	TCCCCGCGAGTCGA GACAC	1032	1206	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	175	TCCTAAGCCCCCAGGT TCTCC	TCCCCGCGAGTCGA GACAC	1032	1206	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	179	TCCTAAGCCCCCAGGT TCTCC	TCCTTCCCCGAGT CGAG	1032	1210	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	179	TCCTAAGCCCCCAGGT TCTCC	TCCTTCCCCGCAGT CGAG	1032	1210	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	179	TCCTAAGCCCCCAGGT TCTCC	TCCTTCCCCGCAGT CGAG	1032	1210	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	179	TCCTAAGCCCCCAGGT TCTCC	TCCTTCCCCGCAGT CGAG	1032	1210	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	179	TCCTAAGCCCCCAGGT TCTCC	TCCTTCCCCGCAGT CGAG	1032	1210	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	179	TCCTAAGCCCCCAGGT TCTCC	TCCTTCCCCGCAGT CGAG	1032	1210	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	183	TCCTAAGCCCCCAGGT TCTCC	TCTCTCCTTCCCCCG CAGTC	1032	1214	no
XMRV gag	83.5	AGCGGCGGACCTC TCATTGACCTT	67	58.3	183	TCCTAAGCCCCCAGGT TCTCC	TCTCTCCTTCCCCCG CAGTC	1032	1214	no
XMRV gag	82.7	TAGCGGCGGACCT CTCATTGACCTT	66.7	56	183	TCCTAAGCCCCCAGGT TCTCC	TCTCTCCTTCCCCCG CAGTC	1032	1214	no
XMRV gag	82.5	AGCGGCGGACCTC TCATTGACCT	66.6	60.9	183	TCCTAAGCCCCCAGGT TCTCC	TCTCTCCTTCCCCCG CAGTC	1032	1214	no
XMRV gag	82.3	ATAGCGGCGGACCT TCTCATTGACCT	66.5	56	183	TCCTAAGCCCCCAGGT TCTCC	TCTCTCCTTCCCCCG CAGTC	1032	1214	no
XMRV gag	81.6	TAGCGGCGGACCT CTCATTGACCT	66.3	58.3	183	TCCTAAGCCCCCAGGT TCTCC	TCTCTCCTTCCCCCG CAGTC	1032	1214	no
XMRV gag	84.9	AGCGGCGGACCTC TCATTGACCTTC	67.6	60	195	TCCTAAGCCCCCAGGT TCTCC	GCTGCGGGGAGGG TCTCTC	1032	1226	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV gag	83.5	AGCGGGGACCTC TCATTGACCTT	67	58.3	195	TCCTAAGCCCCCAGGT TCTCC	GCTGCGGGAGGG TCTCTC	1032	1226	no
XMRV gag	82.7	TAGCGGGGACCT CTCATTGACCTT	66.7	56	195	TCCTAAGCCCCCAGGT TCTCC	GCTGCGGGAGGG TCTCTC	1032	1226	no
XMRV gag	82.5	AGCGGGGACCTC TCATTGACCT	66.6	60.9	195	TCCTAAGCCCCCAGGT TCTCC	GCTGCGGGAGGG TCTCTC	1032	1226	no
XMRV gag	82.3	ATAGCGGGGACCT TCTCATTGACCT	66.5	56	195	TCCTAAGCCCCCAGGT TCTCC	GCTGCGGGAGGG TCTCTC	1032	1226	no
XMRV gag	81.6	TAGCGGGGACCT CTCATTGACCT	66.3	58.3	195	TCCTAAGCCCCCAGGT TCTCC	GCTGCGGGAGGG TCTCTC	1032	1226	no
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	173	CCAGGATCTGAGAG AAGTCAACAAG	AGGCGAAGAGAG GCTGACTG	2919	3091	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	176	CCAGGATCTGAGAG AAGTCAACAAG	CAAAGGCGAAGA GAGGCTGAC	2919	3094	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	122	CAGGATCTGAGAGA AGTCAACAAG	TCCTTTAAATCAAG CACAGTGTACC	2920	3041	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	169	CAGGATCTGAGAGA AGTCAACAAG	CGAAGAGAGGCT GACTGGTG	2920	3088	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	183	CAGGATCTGAGAGA AGTCAACAAG	TCTCCACTCAAAG GCGAAGAG	2920	3102	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	174	AGGATCTGAGAGAA GTCAACAAGC	CAAAGGCGAAGA GAGGCTGAC	2921	3094	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	174	AGGATCTGAGAGAA GTCAACAAGC	CAAAGGCGAAGA GAGGCTGAC	2921	3094	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	181	AGGATCTGAGAGAA GTCAACAAGC	CTCCACTCAAAGG CGAAGAGAG	2921	3101	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	181	AGGATCTGAGAGAA GTCAACAAGC	CTCCACTCAAAGG CGAAGAGAG	2921	3101	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	120	GGATCTGAGAGAAAG TCAACAAGC	TCCTTTAAATCAAG CACAGTGTACC	2922	3041	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	120	GGATCTGAGAGAAAG TCAACAAGC	TCCTTTAAATCAAG CACAGTGTACC	2922	3041	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	123	GGATCTGAGAGAAAG TCAACAAGC	GCATCCTTTAAATC AAGCACAGTG	2922	3044	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	181	GGATCTGAGAGAAAG TCAACAAGC	TCTCCACTCAAAG GCGAAGAG	2922	3102	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	181	GGATCTGAGAGAAAG TCAACAAGC	TCTCCACTCAAAG GCGAAGAG	2922	3102	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	124	GATCTGAGAGAAAGT CAACAAGCG	AGGCATCCTTTAA ATCAAGCACAG	2923	3046	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	128	AAGTCAACAAGCGG GTGGAAG	TCTCAGGCAGAAA AAGGCATCC	2933	3060	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	128	AAGTCAACAAGCGG GTGGAAG	TCTCAGGCAGAAA AAGGCATCC	2933	3060	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	147	AAGTCAACAAGCGG GTGGAAG	GCTGACTGGTGGG GTGGAG	2933	3079	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	147	AAGTCAACAAGCGG GTGGAAG	GCTGACTGGTGGG GTGGAG	2933	3079	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	159	AAGTCAACAAGCGG GTGGAAG	AGGCGAAGAGAG GCTGACTG	2933	3091	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	159	AAGTCAACAAGCGG GTGGAAG	AGGCGAAGAGAG GCTGACTG	2933	3091	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	162	AAGTCAACAAGCGG GTGGAAG	CAAAGGCGAAGA GAGGCTGAC	2933	3094	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	162	AAGTCAACAAGCGG GTGGAAG	CAAAGGCGAAGA GAGGCTGAC	2933	3094	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	108	AGTCAACAAGCGGG TGGAAG	TCCTTTAAATCAAG CACAGTGATACC	2934	3041	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	108	AGTCAACAAGCGGG TGGAAG	TCCTTTAAATCAAG CACAGTGATACC	2934	3041	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	111	AGTCAACAAGCGGG TGGAAG	GCATCCTTTAAATC AAGCACAGTG	2934	3044	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	111	AGTCAACAAGCGGG TGGAAG	GCATCCTTTAAATC AAGCACAGTG	2934	3044	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	113	AGTCAACAAGCGGG TGGAAG	AGGCATCCTTTAA ATCAAGCACAG	2934	3046	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	113	AGTCAACAAGCGGG TGGAAG	AGGCATCCTTTAA ATCAAGCACAG	2934	3046	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	132	AGTCAACAAGCGGG TGGAAG	TGGAGTCTCAGGC AGAAAAAGG	2934	3065	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	132	AGTCAACAAGCGGG TGGAAG	TGGAGTCTCAGGC AGAAAAAGG	2934	3065	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	140	AGTCAACAAGCGGG TGGAAG	TGGTGGGGTGA GTCTCAG	2934	3073	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	140	AGTCAACAAGCGGG TGGAAG	TGGTGGGGTGA GTCTCAG	2934	3073	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	142	AGTCAACAAGCGGG TGGAAG	ACTGGTGGGGTG GAGTCTC	2934	3075	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	142	AGTCAACAAGCGGG TGGAAG	ACTGGTGGGGTG GAGTCTC	2934	3075	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	155	AGTCAACAAGCGGG TGGAAG	CGAAGAGAGGCT GACTGGTG	2934	3088	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	155	AGTCAACAAGCGGG TGGAAG	CGAAGAGAGGCT GACTGGTG	2934	3088	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	200	AGTCAACAAGCGGG TGGAAG	TCAGTTGTCCTGA GATTCCTCATC	2934	3133	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	200	AGTCAACAAGCGGG TGGAAG	TCAGTTGTCCTGA GATTCCTCATC	2934	3133	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	106	TCAACAAGCGGGTG GAAGAC	TCCTTTAAATCAAG CACAGTGATCC	2936	3041	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	106	TCAACAAGCGGGTG GAAGAC	TCCTTTAAATCAAG CACAGTGATCC	2936	3041	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	109	TCAACAAGCGGGTG GAAGAC	GCATCCTTTAAATC AAGCACAGTG	2936	3044	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	109	TCAACAAGCGGGTG GAAGAC	GCATCCTTTAAATC AAGCACAGTG	2936	3044	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	111	TCAACAAGCGGGTG GAAGAC	AGGCATCCTTTAA ATCAAGCACAG	2936	3046	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	111	TCAACAAGCGGGTG GAAGAC	AGGCATCCTTTAA ATCAAGCACAG	2936	3046	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	125	TCAACAAGCGGGTG GAAGAC	TCTCAGGCAGAAA AAGGCATCC	2936	3060	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	125	TCAACAAGCGGGTG GAAGAC	TCTCAGGCAGAAA AAGGCATCC	2936	3060	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	130	TCAACAAGCGGGTG GAAGAC	TGGAGTCTCAGGC AGAAAAAGG	2936	3065	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	138	TCAACAAGCGGGTG GAAGAC	TGGTGGGGTGA GTCTCAG	2936	3073	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	138	TCAACAAGCGGGTG GAAGAC	TGGTGGGGTGA GTCTCAG	2936	3073	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	140	TCAACAAGCGGGTG GAAGAC	ACTGTGGGGTG GAGTCTC	2936	3075	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	140	TCAACAAGCGGGTG GAAGAC	ACTGTGGGGTG GAGTCTC	2936	3075	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	144	TCAACAAGCGGGTG GAAGAC	GCTGACTGGTGGG GTGGAG	2936	3079	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	144	TCAACAAGCGGGTG GAAGAC	GCTGACTGGTGGG GTGGAG	2936	3079	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	153	TCAACAAGCGGGTG GAAGAC	CGAAGAGAGGCT GACTGGTG	2936	3088	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	153	TCAACAAAGCGGGTG GAAGAC	CGAAGAGAGGCT GACTGGTG	2936	3088	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	156	TCAACAAAGCGGGTG GAAGAC	AGGCGAAGAGAG GCTGACTG	2936	3091	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	156	TCAACAAAGCGGGTG GAAGAC	AGGCGAAGAGAG GCTGACTG	2936	3091	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	159	TCAACAAAGCGGGTG GAAGAC	CAAAGGCGAAGA GAGGCTGAC	2936	3094	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	159	TCAACAAAGCGGGTG GAAGAC	CAAAGGCGAAGA GAGGCTGAC	2936	3094	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	167	TCAACAAAGCGGGTG GAAGAC	TCTCCACTCAAAG GCGAAGAG	2936	3102	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	167	TCAACAAAGCGGGTG GAAGAC	TCTCCACTCAAAG GCGAAGAG	2936	3102	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	198	TCAACAAAGCGGGTG GAAGAC	TCAGTTGTCTCTGA GATTCCTCATC	2936	3133	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	198	TCAACAAAGCGGGTG GAAGAC	TCAGTTGTCTCTGA GATTCCTCATC	2936	3133	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	104	AACAAGCGGGTGGA AGACATC	TCCTTTAAATCAAG CACAGTGATCC	2938	3041	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCTTACAACC	68.8	64	104	AACAAGCGGGTGGA AGACATC	TCCTTTAAATCAAG CACAGTGATCC	2938	3041	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	109	AACAAGCGGGTGGA AGACATC	AGGCATCCTTTAA ATCAAGCACAG	2938	3046	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	109	AACAAGCGGGTGGA AGACATC	AGGCATCCTTTAA ATCAAGCACAG	2938	3046	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	123	AACAAGCGGGTGGA AGACATC	TCTCAGGCAGAAA AAGGCATCC	2938	3060	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	123	AACAAGCGGGTGGA AGACATC	TCTCAGGCAGAAA AAGGCATCC	2938	3060	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	128	AACAAGCGGGTGGA AGACATC	TGGAGTCTCAGGC AGAAAAAGG	2938	3065	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	142	AACAAGCGGGTGGA AGACATC	GCTGACTGGTGGG GTGGAG	2938	3079	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	142	AACAAGCGGGTGGA AGACATC	GCTGACTGGTGGG GTGGAG	2938	3079	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	154	AACAAGCGGGTGGA AGACATC	AGGCGAAGAGAG GCTGACTG	2938	3091	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	154	AACAAGCGGGTGGA AGACATC	AGGCGAAGAGAG GCTGACTG	2938	3091	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	157	AACAAGCGGGTGGA AGACATC	CAAGGCGAAGA GAGGCTGAC	2938	3094	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	157	AACAAGCGGGTGGA AGACATC	CAAGGCGAAGA GAGGCTGAC	2938	3094	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	184	AACAAGCGGGTGGA AGACATC	AGATTCCTCTCT GGATCTCTCC	2938	3121	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	184	AACAAGCGGGTGGA AGACATC	AGATTCCTCTCT GGATCTCTCC	2938	3121	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	106	ACAAGCGGGTGGAA GACATC	GCATCCTTTAAATC AAGCACAGTG	2939	3044	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	106	ACAAGCGGGTGGAA GACATC	GCATCCTTTAAATC AAGCACAGTG	2939	3044	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	135	ACAAGCGGGTGGAA GACATC	TGGTGGGGTGGAA GTCTCAG	2939	3073	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	135	ACAAGCGGGTGGAA GACATC	TGGTGGGGTGGAA GTCTCAG	2939	3073	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	137	ACAAGCGGGTGGAA GACATC	ACTGGTGGGGTG GAGTCTC	2939	3075	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	137	ACAAGCGGGTGGAA GACATC	ACTGGTGGGGTG GAGTCTC	2939	3075	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	150	ACAAGCGGGTGGAA GACATC	CGAAGAGAGGCT GACTGGTG	2939	3088	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	150	ACAAGCGGGTGGAA GACATC	CGAAGAGAGGCT GACTGGTG	2939	3088	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	164	ACAAGCGGGTGGAA GACATC	TCTCCACTCAAAG GCGAAGAG	2939	3102	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	164	ACAAGCGGGTGGAA GACATC	TCTCCACTCAAAG GCGAAGAG	2939	3102	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	195	ACAAGCGGGTGGAA GACATC	TCAGTTGTCTCTGA GATTCCCATC	2939	3133	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	195	ACAAGCGGGTGGAA GACATC	TCAGTTGTCTCTGA GATTCCCATC	2939	3133	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	198	ACAAGCGGGTGGAA GACATC	AGGTCAGTTGTCC TGAGATTCC	2939	3136	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	198	ACAAGCGGGTGGAA GACATC	AGGTCAGTTGTCC TGAGATTCC	2939	3136	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	102	CAAGCGGGTGGAAAG ACATCC	TCCTTTAAATCAAG CACAGTGATCC	2940	3041	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	102	CAAGCGGGTGGAAAG ACATCC	TCCTTTAAATCAAG CACAGTGATCC	2940	3041	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	107	CAAGCGGGTGGAAAG ACATCC	AGGCATCCTTTAA ATCAAGCACAG	2940	3046	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	107	CAAGCGGGTGGAAAG ACATCC	AGGCATCCTTTAA ATCAAGCACAG	2940	3046	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	152	CAAGCGGGTGGAAAG ACATCC	AGGCGAAGAGAG GCTGACTG	2940	3091	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	152	CAAGCGGGTGGAAAG ACATCC	AGGCGAAGAGAG GCTGACTG	2940	3091	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	155	CAAGCGGGTGGAAAG ACATCC	CAAAGGCGAAGA GAGGCTGAC	2940	3094	yes
XMRV pol	75.1	CCCACCGTGCCCA ACCCTTACAACC	68.8	64	155	CAAGCGGGTGGAAAG ACATCC	CAAAGGCGAAGA GAGGCTGAC	2940	3094	yes
XMRV pol	75.6	ACCGTGCCCAACC CTTACAACCTCT	67.1	56	148	AAGCGGGTGGAAAG CATCC	CGAAGAGAGGCT GACTGGTG	2941	3088	yes
XMRV pol	77.2	CCCACCAAGTCAGC CTCTCTTCGCCT	68.9	64	116	TGATTTAAAGGATGC CTTTTCTGC	GTCTGGTCCAGGT CAGTTGTC	3030	3145	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	116	TGATTTAAAGGATGC CTTTTCTGC	GTCTGGTCCAGGT CAGTTGTC	3030	3145	yes
XMRV pol	76	CCACCAGTCAGC TCTCTTCGCCT	67.3	60	116	TGATTTAAAGGATGC CTTTTCTGC	GTCTGGTCCAGGT CAGTTGTC	3030	3145	yes
XMRV pol	75.5	CCACCAGTCAGC TCTCTTCGCCT	67	62.5	116	TGATTTAAAGGATGC CTTTTCTGC	GTCTGGTCCAGGT CAGTTGTC	3030	3145	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	126	AAGGATGCCTTTTCT TGCCTGAG	TTTGAAACCCTGT GGGAGTCTG	3037	3162	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	165	AAGGATGCCTTTTCT TGCCTGAG	GTCTCTGTGCAGT GCCTCATC	3037	3201	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	84	AGGATGCCTTTTCT GCCTGAG	AGATTCCTCATCTCT GGATCTCTCC	3038	3121	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	84	AGGATGCCTTTTCT GCCTGAG	AGATTCCTCATCTCT GGATCTCTCC	3038	3121	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	96	AGGATGCCTTTTCT GCCTGAG	TCAGTTGTCCTGA GATTCCTCATC	3038	3133	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	96	AGGATGCCTTTTCT GCCTGAG	TCAGTTGTCCTGA GATTCCTCATC	3038	3133	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	126	AAGGATGCCTTTTCT TGCCTGAG	TTTGAAACCCTGT GGGAGTCTG	3038	3162	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	149	AGGATGCCTTTTCT GCCTGAG	CTCATCAAAACAGG GTGGGACTG	3038	3186	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	149	AGGATGCCTTTTCT GCCTGAG	CTCATCAAAACAGG GTGGGACTG	3038	3186	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCT	67.3	60	149	AGGATGCCTTTTCT GCCTGAG	CTCATCAAAACAGG GTGGGACTG	3038	3186	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	149	AGGATGCCTTTTCT GCCTGAG	CTCATCAAAACAGG GTGGGACTG	3038	3186	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	165	AAGGATGCCTTTTCT TGCCTGAG	GTCTCTGTGCAGT GCCTCATC	3038	3201	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	178	AGGATGCCTTTTCT GCCTGAG	CGGAAATCTGCTA GGTCTCTGTG	3038	3215	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCCT	67.7	66.7	178	AGGATGCCTTTTCT GCCTGAG	CGGAAATCTGCTA GGTCTCTGTG	3038	3215	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCT	67.3	60	178	AGGATGCCTTTTCT GCCTGAG	CGGAAATCTGCTA GGTCTCTGTG	3038	3215	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	178	AGGATGCCTTTTCT GCCTGAG	CGGAAATCTGCTA GGTCTCTGTG	3038	3215	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	189	AGGATGCCTTTTCT GCCTGAG	GGTGTGGATCCG GAAATCTG	3038	3226	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCCT	67.7	66.7	189	AGGATGCCTTTTCT GCCTGAG	GGTGTGGATCCG GAAATCTG	3038	3226	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	98	GGATGCCTTTTCTG CCTGAG	AGGTCAGTTGTCC TGAGATTCC	3039	3136	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCCT	67.7	66.7	98	GGATGCCTTTTCTG CCTGAG	AGGTCAGTTGTCC TGAGATTCC	3039	3136	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCT	67.3	60	98	GGATGCCTTTTCTG CCTGAG	AGGTCAGTTGTCC TGAGATTCC	3039	3136	yes
XMRV pol	75.5	CCACCAGTCAGCC	67	62.5	98	GGATGCCTTTTCTG	AGGTCAGTTGTCC	3039	3136	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
		TCTCTTCGCCT				CCTGAG	TGAGATTCC			
XMRV pol	77.2	CCCACCAAGTCAGC CTCTCTTCGCCT	68.9	64	107	GGATGCCCTTTTCTG CCTGAG	GTCTGGTCCAGGT CAGTTGTC	3039	3145	yes
XMRV pol	76.5	CCCACCAAGTCAGC CTCTCTTCGCC	67.7	66.7	107	GGATGCCCTTTTCTG CCTGAG	GTCTGGTCCAGGT CAGTTGTC	3039	3145	yes
XMRV pol	76	CCACCAAGTCAGCC TCTCTTCGCCT	67.3	60	107	GGATGCCCTTTTCTG CCTGAG	GTCTGGTCCAGGT CAGTTGTC	3039	3145	yes
XMRV pol	75.5	CCACCAAGTCAGCC TCTCTTCGCCT	67	62.5	107	GGATGCCCTTTTCTG CCTGAG	GTCTGGTCCAGGT CAGTTGTC	3039	3145	yes
XMRV pol	77.2	CCCACCAAGTCAGC CTCTCTTCGCCT	68.9	64	109	GGATGCCCTTTTCTG CCTGAG	GAGTCTGGTCCAG GTCAGTTG	3039	3147	no
XMRV pol	76.5	CCCACCAAGTCAGC CTCTCTTCGCC	67.7	66.7	109	GGATGCCCTTTTCTG CCTGAG	GAGTCTGGTCCAG GTCAGTTG	3039	3147	no
XMRV pol	76	CCACCAAGTCAGCC TCTCTTCGCCT	67.3	60	109	GGATGCCCTTTTCTG CCTGAG	GAGTCTGGTCCAG GTCAGTTG	3039	3147	no
XMRV pol	75.5	CCACCAAGTCAGCC TCTCTTCGCCT	67	62.5	109	GGATGCCCTTTTCTG CCTGAG	GAGTCTGGTCCAG GTCAGTTG	3039	3147	no
XMRV pol	77.2	CCCACCAAGTCAGC CTCTCTTCGCCT	68.9	64	119	GGATGCCCTTTTCTG CCTGAG	AACCTGTGGGAG TCTGGTC	3039	3157	yes
XMRV pol	76.5	CCCACCAAGTCAGC CTCTCTTCGCC	67.7	66.7	119	GGATGCCCTTTTCTG CCTGAG	AACCTGTGGGAG TCTGGTC	3039	3157	yes
XMRV pol	76	CCACCAAGTCAGCC TCTCTTCGCCT	67.3	60	119	GGATGCCCTTTTCTG CCTGAG	AACCTGTGGGAG TCTGGTC	3039	3157	yes
XMRV pol	75.5	CCACCAAGTCAGCC TCTCTTCGCCT	67	62.5	119	GGATGCCCTTTTCTG CCTGAG	AACCTGTGGGAG TCTGGTC	3039	3157	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	124	GGATGCCCTTTTCTGCTGAG	TTTGAAACCCCTGTGGGAGTCTG	3039	3162	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	124	GGATGCCCTTTTCTGCTGAG	TTTGAAACCCCTGTGGGAGTCTG	3039	3162	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	124	GGATGCCCTTTTCTGCTGAG	TTTGAAACCCCTGTGGGAGTCTG	3039	3162	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	124	GGATGCCCTTTTCTGCTGAG	TTTGAAACCCCTGTGGGAGTCTG	3039	3162	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	157	GGATGCCCTTTTCTGCTGAG	GTGCAGTGCCTCATCAACACAG	3039	3195	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	157	GGATGCCCTTTTCTGCTGAG	GTGCAGTGCCTCATCAACACAG	3039	3195	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	157	GGATGCCCTTTTCTGCTGAG	GTGCAGTGCCTCATCAACACAG	3039	3195	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	157	GGATGCCCTTTTCTGCTGAG	GTGCAGTGCCTCATCAACACAG	3039	3195	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	159	GGATGCCCTTTTCTGCTGAG	CTGTGCAGTGCCTCATCAAC	3039	3197	yes
XMRV pol	76.5	CCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	159	GGATGCCCTTTTCTGCTGAG	CTGTGCAGTGCCTCATCAAC	3039	3197	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	159	GGATGCCCTTTTCTGCTGAG	CTGTGCAGTGCCTCATCAAC	3039	3197	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	159	GGATGCCCTTTTCTGCTGAG	CTGTGCAGTGCCTCATCAAC	3039	3197	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	163	GGATGCCCTTTTCTGCTGAG	GTCTCTGTGCAGTGCCTCATC	3039	3201	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	163	GGATGCCCTTTTCTGCTGAG	GTCTCTGTGCAGTGCCTCATC	3039	3201	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCT	67.3	60	163	GGATGCCCTTTTCTGCTGAG	GTCTCTGTGCAGTGCCTCATC	3039	3201	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	163	GGATGCCCTTTTCTGCTGAG	GTCTCTGTGCAGTGCCTCATC	3039	3201	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	168	GGATGCCCTTTTCTGCTGAG	GCTAGGTCTCTGTGCAGTGC	3039	3206	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	168	GGATGCCCTTTTCTGCTGAG	GCTAGGTCTCTGTGCAGTGC	3039	3206	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCT	67.3	60	168	GGATGCCCTTTTCTGCTGAG	GCTAGGTCTCTGTGCAGTGC	3039	3206	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	168	GGATGCCCTTTTCTGCTGAG	GCTAGGTCTCTGTGCAGTGC	3039	3206	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	178	GGATGCCCTTTTCTGCTGAG	CCGGAAATCTGCTAGGTCTCTG	3039	3216	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	178	GGATGCCCTTTTCTGCTGAG	CCGGAAATCTGCTAGGTCTCTG	3039	3216	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCT	67.3	60	178	GGATGCCCTTTTCTGCTGAG	CCGGAAATCTGCTAGGTCTCTG	3039	3216	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	178	GGATGCCCTTTTCTGCTGAG	CCGGAAATCTGCTAGGTCTCTG	3039	3216	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	181	GGATGCCTTTTCTG CCTGAG	GATCCGGAAATCT GCTAGGTCTC	3039	3219	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	181	GGATGCCTTTTCTG CCTGAG	GATCCGGAAATCT GCTAGGTCTC	3039	3219	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCT	67.3	60	181	GGATGCCTTTTCTG CCTGAG	GATCCGGAAATCT GCTAGGTCTC	3039	3219	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	181	GGATGCCTTTTCTG CCTGAG	GATCCGGAAATCT GCTAGGTCTC	3039	3219	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	198	GGATGCCTTTTCTG CCTGAG	ATCAAGTCTGGGT GCTGGATC	3039	3236	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	198	GGATGCCTTTTCTG CCTGAG	ATCAAGTCTGGGT GCTGGATC	3039	3236	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCT	67.3	60	198	GGATGCCTTTTCTG CCTGAG	ATCAAGTCTGGGT GCTGGATC	3039	3236	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	198	GGATGCCTTTTCTG CCTGAG	ATCAAGTCTGGGT GCTGGATC	3039	3236	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	82	GATGCCTTTTCTGC CTGAGAC	AGATCCCACATCTCT GGATCTCTCC	3040	3121	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	82	GATGCCTTTTCTGC CTGAGAC	AGATCCCACATCTCT GGATCTCTCC	3040	3121	no
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	94	GATGCCTTTTCTGC CTGAGAC	TCAGTTGTCCTGA GATCCCACATC	3040	3133	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	94	GATGCCTTTTCTGC CTGAGAC	TCAGTTGTCCTGA GATCCCACATC	3040	3133	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCT	67.3	60	94	GATGCCTTTTCTGC CTGAGAC	TCAGTTGTCTGA GATCCCATC	3040	3133	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	94	GATGCCTTTTCTGC CTGAGAC	TCAGTTGTCTGA GATCCCATC	3040	3133	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	110	GATGCCTTTTCTGC CTGAGAC	GGGAGTCTGGTCC AGGTCAG	3040	3149	no
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	110	GATGCCTTTTCTGC CTGAGAC	GGGAGTCTGGTCC AGGTCAG	3040	3149	no
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCT	67.3	60	110	GATGCCTTTTCTGC CTGAGAC	GGGAGTCTGGTCC AGGTCAG	3040	3149	no
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	110	GATGCCTTTTCTGC CTGAGAC	GGGAGTCTGGTCC AGGTCAG	3040	3149	no
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	126	GATGCCTTTTCTGC CTGAGAC	GTTTTGAAACCCT GTGGGAGTC	3040	3165	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	126	GATGCCTTTTCTGC CTGAGAC	GTTTTGAAACCCT GTGGGAGTC	3040	3165	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCT	67.3	60	126	GATGCCTTTTCTGC CTGAGAC	GTTTTGAAACCCT GTGGGAGTC	3040	3165	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	126	GATGCCTTTTCTGC CTGAGAC	GTTTTGAAACCCT GTGGGAGTC	3040	3165	yes
XMRV pol	77.2	CCCACCAGTCAGC CTCTCTTCGCCT	68.9	64	132	GATGCCTTTTCTGC CTGAGAC	GGGACTGTTTTG AAACCCGTG	3040	3171	yes
XMRV pol	76.5	CCCACCAGTCAGC CTCTCTTCGCC	67.7	66.7	132	GATGCCTTTTCTGC CTGAGAC	GGGACTGTTTTG AAACCCGTG	3040	3171	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	147	GATGCCCTTTTCTGCTGAGAC	CTCATCAAAACAGG GTGGGACTG	3040	3186	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	147	GATGCCCTTTTCTGCTGAGAC	CTCATCAAAACAGG GTGGGACTG	3040	3186	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	147	GATGCCCTTTTCTGCTGAGAC	CTCATCAAAACAGG GTGGGACTG	3040	3186	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	147	GATGCCCTTTTCTGCTGAGAC	CTCATCAAAACAGG GTGGGACTG	3040	3186	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	181	GATGCCCTTTTCTGCTGAGAC	GGATCCGGAAATC TGCTAGGTC	3040	3220	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	181	GATGCCCTTTTCTGCTGAGAC	GGATCCGGAAATC TGCTAGGTC	3040	3220	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	181	GATGCCCTTTTCTGCTGAGAC	GGATCCGGAAATC TGCTAGGTC	3040	3220	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	181	GATGCCCTTTTCTGCTGAGAC	GGATCCGGAAATC TGCTAGGTC	3040	3220	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	183	GATGCCCTTTTCTGCTGAGAC	CTGGATCCGGAAA TCTGCTAGG	3040	3222	yes
XMRV pol	76.5	CCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	183	GATGCCCTTTTCTGCTGAGAC	CTGGATCCGGAAA TCTGCTAGG	3040	3222	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	183	GATGCCCTTTTCTGCTGAGAC	CTGGATCCGGAAA TCTGCTAGG	3040	3222	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	183	GATGCCCTTTTCTGCTGAGAC	CTGGATCCGGAAA TCTGCTAGG	3040	3222	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	187	GATGCCCTTTTCTGCTGAGAC	GGTGCTGGATCCG GAAATCTG	3040	3226	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	187	GATGCCCTTTTCTGCTGAGAC	GGTGCTGGATCCG GAAATCTG	3040	3226	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	187	GATGCCCTTTTCTGCTGAGAC	GGTGCTGGATCCG GAAATCTG	3040	3226	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	187	GATGCCCTTTTCTGCTGAGAC	GGTGCTGGATCCG GAAATCTG	3040	3226	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	188	GATGCCCTTTTCTGCTGAGAC	GGTGCTGGATCC GGAAATC	3040	3227	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	188	GATGCCCTTTTCTGCTGAGAC	GGTGCTGGATCC GGAAATC	3040	3227	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	188	GATGCCCTTTTCTGCTGAGAC	GGTGCTGGATCC GGAAATC	3040	3227	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	188	GATGCCCTTTTCTGCTGAGAC	GGTGCTGGATCC GGAAATC	3040	3227	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	200	GATGCCCTTTTCTGCTGAGAC	AGGATCAAGTCTG GGTGCTG	3040	3239	yes
XMRV pol	76.5	CCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	200	GATGCCCTTTTCTGCTGAGAC	AGGATCAAGTCTG GGTGCTG	3040	3239	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	200	GATGCCCTTTTCTGCTGAGAC	AGGATCAAGTCTG GGTGCTG	3040	3239	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	200	GATGCCCTTTTCTGCTGAGAC	AGGATCAAGTCTG GGTGCTG	3040	3239	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	96	ATGCCCTTTTCTGCCTGAGAC	AGGTCAGTTGTCCTGAGATTCC	3041	3136	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	96	ATGCCCTTTTCTGCCTGAGAC	AGGTCAGTTGTCCTGAGATTCC	3041	3136	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	96	ATGCCCTTTTCTGCCTGAGAC	AGGTCAGTTGTCCTGAGATTCC	3041	3136	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	96	ATGCCCTTTTCTGCCTGAGAC	AGGTCAGTTGTCCTGAGATTCC	3041	3136	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	107	ATGCCCTTTTCTGCCTGAGAC	GAGTCTGGTCCAGGTCAGTTG	3041	3147	no
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	107	ATGCCCTTTTCTGCCTGAGAC	GAGTCTGGTCCAGGTCAGTTG	3041	3147	no
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	107	ATGCCCTTTTCTGCCTGAGAC	GAGTCTGGTCCAGGTCAGTTG	3041	3147	no
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	107	ATGCCCTTTTCTGCCTGAGAC	GAGTCTGGTCCAGGTCAGTTG	3041	3147	no
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	196	ATGCCCTTTTCTGCCTGAGAC	ATCAAGTCTGGTGCTGGATC	3041	3236	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	196	ATGCCCTTTTCTGCCTGAGAC	ATCAAGTCTGGTGCTGGATC	3041	3236	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	196	ATGCCCTTTTCTGCCTGAGAC	ATCAAGTCTGGTGCTGGATC	3041	3236	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	196	ATGCCCTTTTCTGCCTGAGAC	ATCAAGTCTGGTGCTGGATC	3041	3236	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	80	TGCCTTTTCTGCCTGAGACTC	AGATTCCCATTCTCTGGATCTCTCC	3042	3121	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	80	TGCCTTTTCTGCCTGAGACTC	AGATTCCCATTCTCTGGATCTCTCC	3042	3121	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	80	TGCCTTTTCTGCCTGAGACTC	AGATTCCCATTCTCTGGATCTCTCC	3042	3121	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	80	TGCCTTTTCTGCCTGAGACTC	AGATTCCCATTCTCTGGATCTCTCC	3042	3121	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	92	TGCCTTTTCTGCCTGAGACTC	TCAGTTGTCTCTGAGATTCCCATC	3042	3133	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	92	TGCCTTTTCTGCCTGAGACTC	TCAGTTGTCTCTGAGATTCCCATC	3042	3133	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	92	TGCCTTTTCTGCCTGAGACTC	TCAGTTGTCTCTGAGATTCCCATC	3042	3133	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	92	TGCCTTTTCTGCCTGAGACTC	TCAGTTGTCTCTGAGATTCCCATC	3042	3133	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	104	TGCCTTTTCTGCCTGAGACTC	GTCTGGTCCAGGTGAGTTGTC	3042	3145	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	104	TGCCTTTTCTGCCTGAGACTC	GTCTGGTCCAGGTGAGTTGTC	3042	3145	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	104	TGCCTTTTCTGCCTGAGACTC	GTCTGGTCCAGGTGAGTTGTC	3042	3145	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	104	TGCCTTTTCTGCCTGAGACTC	GTCTGGTCCAGGTGAGTTGTC	3042	3145	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	130	TGCCTTTTCTGCCTGAGACTC	GGGACTGTTTTTGAAACCCCTGTG	3042	3171	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	130	TGCCTTTTCTGCCTGAGACTC	GGGACTGTTTTTGAAACCCCTGTG	3042	3171	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	130	TGCCTTTTCTGCCTGAGACTC	GGGACTGTTTTTGAAACCCCTGTG	3042	3171	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	130	TGCCTTTTCTGCCTGAGACTC	GGGACTGTTTTTGAAACCCCTGTG	3042	3171	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	145	TGCCTTTTCTGCCTGAGACTC	CTCATCAAAACAGG GTGGGACTG	3042	3186	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	145	TGCCTTTTCTGCCTGAGACTC	CTCATCAAAACAGG GTGGGACTG	3042	3186	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	145	TGCCTTTTCTGCCTGAGACTC	CTCATCAAAACAGG GTGGGACTG	3042	3186	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	145	TGCCTTTTCTGCCTGAGACTC	CTCATCAAAACAGG GTGGGACTG	3042	3186	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	154	TGCCTTTTCTGCCTGAGACTC	GTGCAGTGCCTCA TCAAAACAG	3042	3195	yes
XMRV pol	76.5	CCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	154	TGCCTTTTCTGCCTGAGACTC	GTGCAGTGCCTCA TCAAAACAG	3042	3195	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	154	TGCCTTTTCTGCCTGAGACTC	GTGCAGTGCCTCA TCAAAACAG	3042	3195	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	154	TGCCTTTTCTGCCTGAGACTC	GTGCAGTGCCTCA TCAAAACAG	3042	3195	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	156	TGCCCTTTTCTGCCTGAGACTC	CTGTGCAGTGCCTCATCAAAAC	3042	3197	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	156	TGCCCTTTTCTGCCTGAGACTC	CTGTGCAGTGCCTCATCAAAAC	3042	3197	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	156	TGCCCTTTTCTGCCTGAGACTC	CTGTGCAGTGCCTCATCAAAAC	3042	3197	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	156	TGCCCTTTTCTGCCTGAGACTC	CTGTGCAGTGCCTCATCAAAAC	3042	3197	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	160	TGCCCTTTTCTGCCTGAGACTC	GTCTCTGTGCAGTGCCTCATC	3042	3201	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	160	TGCCCTTTTCTGCCTGAGACTC	GTCTCTGTGCAGTGCCTCATC	3042	3201	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	160	TGCCCTTTTCTGCCTGAGACTC	GTCTCTGTGCAGTGCCTCATC	3042	3201	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	160	TGCCCTTTTCTGCCTGAGACTC	GTCTCTGTGCAGTGCCTCATC	3042	3201	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	174	TGCCCTTTTCTGCCTGAGACTC	CGGAAATCTGCTAGGTCTCTGTG	3042	3215	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	174	TGCCCTTTTCTGCCTGAGACTC	CGGAAATCTGCTAGGTCTCTGTG	3042	3215	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	174	TGCCCTTTTCTGCCTGAGACTC	CGGAAATCTGCTAGGTCTCTGTG	3042	3215	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	174	TGCCCTTTTCTGCCTGAGACTC	CGGAAATCTGCTAGGTCTCTGTG	3042	3215	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	175	TGCCCTTTTCTGCCTGAGACTC	CCGGAAATCTGCTAGGTCTCTG	3042	3216	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	175	TGCCCTTTTCTGCCTGAGACTC	CCGGAAATCTGCTAGGTCTCTG	3042	3216	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	175	TGCCCTTTTCTGCCTGAGACTC	CCGGAAATCTGCTAGGTCTCTG	3042	3216	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	175	TGCCCTTTTCTGCCTGAGACTC	CCGGAAATCTGCTAGGTCTCTG	3042	3216	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	179	TGCCCTTTTCTGCCTGAGACTC	GGATCCGGAAATCTGCTAGGTC	3042	3220	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	179	TGCCCTTTTCTGCCTGAGACTC	GGATCCGGAAATCTGCTAGGTC	3042	3220	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	179	TGCCCTTTTCTGCCTGAGACTC	GGATCCGGAAATCTGCTAGGTC	3042	3220	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	179	TGCCCTTTTCTGCCTGAGACTC	GGATCCGGAAATCTGCTAGGTC	3042	3220	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	181	TGCCCTTTTCTGCCTGAGACTC	CTGGATCCGGAAATCTGCTAGG	3042	3222	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	181	TGCCCTTTTCTGCCTGAGACTC	CTGGATCCGGAAATCTGCTAGG	3042	3222	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	181	TGCCCTTTTCTGCCTGAGACTC	CTGGATCCGGAAATCTGCTAGG	3042	3222	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	181	TGCCCTTTTCTGCCTGAGACTC	CTGGATCCGGAAATCTGCTAGG	3042	3222	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	185	TGCCCTTTTCTGCCTGAGACTC	GGTGCTGGATCCG GAAATCTG	3042	3226	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	185	TGCCCTTTTCTGCCTGAGACTC	GGTGCTGGATCCG GAAATCTG	3042	3226	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	185	TGCCCTTTTCTGCCTGAGACTC	GGTGCTGGATCCG GAAATCTG	3042	3226	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	185	TGCCCTTTTCTGCCTGAGACTC	GGTGCTGGATCCG GAAATCTG	3042	3226	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	186	TGCCCTTTTCTGCCTGAGACTC	GGTGCTGGATCC GGAATC	3042	3227	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	186	TGCCCTTTTCTGCCTGAGACTC	GGTGCTGGATCC GGAATC	3042	3227	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	186	TGCCCTTTTCTGCCTGAGACTC	GGTGCTGGATCC GGAATC	3042	3227	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	186	TGCCCTTTTCTGCCTGAGACTC	GGTGCTGGATCC GGAATC	3042	3227	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	94	GCCTTTTCTGCCTGAGACTC	AGGTCAGTTGTCC TGAGATTCC	3043	3136	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	94	GCCTTTTCTGCCTGAGACTC	AGGTCAGTTGTCC TGAGATTCC	3043	3136	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	94	GCCTTTTCTGCCTGAGACTC	AGGTCAGTTGTCC TGAGATTCC	3043	3136	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	94	GCCTTTTCTGCCTGAGACTC	AGGTCAGTTGTCC TGAGATTCC	3043	3136	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	194	GCCTTTTCTGCCTGAGACTC	ATCAAGTCTGGGTGCTGGATC	3043	3236	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	194	GCCTTTTCTGCCTGAGACTC	ATCAAGTCTGGGTGCTGGATC	3043	3236	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	194	GCCTTTTCTGCCTGAGACTC	ATCAAGTCTGGGTGCTGGATC	3043	3236	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	194	GCCTTTTCTGCCTGAGACTC	ATCAAGTCTGGGTGCTGGATC	3043	3236	yes
XMRV pol	77.2	CCCACCAGTCAGCCTCTCTTCGCCT	68.9	64	197	GCCTTTTCTGCCTGAGACTC	AGGATCAAGTCTGGTGCTG	3043	3239	yes
XMRV pol	76.5	CCCACCAGTCAGCCTCTCTTCGCC	67.7	66.7	197	GCCTTTTCTGCCTGAGACTC	AGGATCAAGTCTGGTGCTG	3043	3239	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	89	CTTTTCTGCCTGAGACTCCAC	TCAGTTGTCCTGAGATTCCCATC	3045	3133	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	89	CTTTTCTGCCTGAGACTCCAC	TCAGTTGTCCTGAGATTCCCATC	3045	3133	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	92	CTTTTCTGCCTGAGACTCCAC	AGGTCAGTTGTCC TGAGATTCC	3045	3136	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	92	CTTTTCTGCCTGAGACTCCAC	AGGTCAGTTGTCC TGAGATTCC	3045	3136	yes
XMRV pol	76	CCACCAGTCAGCCTCTCTTCGCCTT	67.3	60	101	CTTTTCTGCCTGAGACTCCAC	GTCTGGTCCAGGT CAGTTGTC	3045	3145	yes
XMRV pol	75.5	CCACCAGTCAGCCTCTCTTCGCCT	67	62.5	101	CTTTTCTGCCTGAGACTCCAC	GTCTGGTCCAGGT CAGTTGTC	3045	3145	yes

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCTT	67.3	60	151	CTTTTCTGCCTGAG ACTCCAC	GTGCAGTGCCTCA TCAAAACAG	3045	3195	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	151	CTTTTCTGCCTGAG ACTCCAC	GTGCAGTGCCTCA TCAAAACAG	3045	3195	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCTT	67.3	60	153	CTTTTCTGCCTGAG ACTCCAC	CTGTGCAGTGCCT CATCAAAAC	3045	3197	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	153	CTTTTCTGCCTGAG ACTCCAC	CTGTGCAGTGCCT CATCAAAAC	3045	3197	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCTT	67.3	60	158	CTTTTCTGCCTGA GACTCCAC	GTCTCTGTGCAGT GCCTCATC	3045	3201	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	158	CTTTTCTGCCTGA GACTCCAC	GTCTCTGTGCAGT GCCTCATC	3045	3201	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCTT	67.3	60	172	CTTTTCTGCCTGA GACTCCAC	CGGAAATCTGCTA GGTCTCTGTG	3045	3215	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	172	CTTTTCTGCCTGA GACTCCAC	CGGAAATCTGCTA GGTCTCTGTG	3045	3215	yes
XMRV pol	76	CCACCAGTCAGCC TCTCTTCGCCTT	67.3	60	192	CTTTTCTGCCTGAG ACTCCAC	ATCAAGTCTGGGT GCTGGATC	3045	3236	yes
XMRV pol	75.5	CCACCAGTCAGCC TCTCTTCGCCT	67	62.5	192	CTTTTCTGCCTGAG ACTCCAC	ATCAAGTCTGGGT GCTGGATC	3045	3236	yes
XRMV env	76.9	TGACCTGTTCTCT CTGACCCGCCA	68.1	60	94	AAACATGGGGACTA AGACTGTATCG	ACGGGGGCCCTA CATTG	6367	6460	no
XRMV env	76.9	TGACCTGTTCTCT CTGACCCGCCA	68.1	60	135	AACATGGGGACTAA GACTGTATCG	GGGGGTAGTGTT CAGTGATC	6368	6502	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	172	AACATGGGGACTAA GACTGTATCG	AGGAGTCCTGGG GAGCATG	6368	6539	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	200	AACATGGGGACTAA GACTGTATCG	TAGAGCCCGCGCC TGAAG	6368	6567	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	90	ACATGGGGACTAAG ACTGTATCG	GCGGGGCCCTACA TTGAG	6369	6458	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	115	ACATGGGGACTAAG ACTGTATCG	TCACGGGATTAGG CCCAATG	6369	6483	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	122	ACATGGGGACTAAG ACTGTATCG	TCAGTGATCACGG GATTAGGC	6369	6490	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	126	ACATGGGGACTAAG ACTGTATCG	CTGTTCAGTGATC ACGGGATTAG	6369	6494	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	155	ACATGGGGACTAAG ACTGTATCG	ATGATCTGCACGG GTTGGG	6369	6523	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	159	AACATGGGGACTAA GACTGTATCG	AGCATGATCTGCA CGGGTTG	6369	6526	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	167	ACATGGGGACTAAG ACTGTATCG	GTCCTGGGGAGCA TGATCTG	6369	6535	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	169	ACATGGGGACTAAG ACTGTATCG	GAGTCCTGGGGA GCATGATC	6369	6537	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	189	ACATGGGGACTAAG ACTGTATCG	GCCTGAAGGAGG AGGACGAG	6369	6557	no
XRMV env	76.9	TGACCTGTCTCTCTGACCCGCCA	68.1	60	195	AAACATGGGGACTA AGACTGTATCG	CCGCGCCTGAAGG AGGAG	6369	6561	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	193	ACATGGGGACTAAGACTGTATCG	CCGCGCCTGAAGGAGGAG	6369	6561	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	189	CATGGGGACTAAGACTGTATCG	CGCCTGAAGGAGGAGGAG	6370	6558	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	153	ATGGGGACTAAGACTGTATCGATC	ATGATCTGCACGGGTTGGG	6371	6523	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	182	ATGGGGACTAAGACTGTATCGATC	AAGGAGGAGGACGAGGAGTC	6371	6552	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	180	GGGGACTAAGACTGTATCGATCC	AAGGAGGAGGACGAGGAGTC	6373	6552	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	86	GGACTAAGACTGTATCGATCCACTG	ACGCGGGGCCCTACATTG	6375	6460	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	109	GGACTAAGACTGTATCGATCCACTG	TCACGGGATTAGGCCCAATG	6375	6483	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	149	GGACTAAGACTGTATCGATCCACTG	ATGATCTGCACGGGTTGGG	6375	6523	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	152	GGACTAAGACTGTATCGATCCACTG	AGCATGATCTGCACGGGTTG	6375	6526	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	178	GGACTAAGACTGTATCGATCCACTG	AAGGAGGAGGACGAGGAGTC	6375	6552	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	187	GGACTAAGACTGTATCGATCCACTG	CCGCGCCTGAAGGAGGAG	6375	6561	no
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	193	GGACTAAGACTGTATCGATCCACTG	TAGAGGCCGCGCCGGAAG	6375	6567	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	83	GACTAAGACTGTATC GATCCACTG	GCGGGCCCTACATTGAG	6376	6458	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	108	GACTAAGACTGTATC GATCCACTG	TCACGGGATTAGGCCCAATG	6376	6483	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	148	GACTAAGACTGTATC GATCCACTG	ATGATCTGCACGGGTTGGG	6376	6523	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	162	GACTAAGACTGTATC GATCCACTG	GAGTCCTGGGGAGCATGATC	6376	6537	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	183	GACTAAGACTGTATC GATCCACTG	CGCTGAAGGAGGAGGAG	6376	6558	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	150	ACTAAGACTGTATCG ATCCACTGG	AGCATGATCTGCA CGGGTTG	6377	6526	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	159	ACTAAGACTGTATCG ATCCACTGG	GTCCTGGGGAGCATGATCTG	6377	6535	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	163	ACTAAGACTGTATCG ATCCACTGG	AGGAGTCCTGGGAGCATG	6377	6539	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	176	ACTAAGACTGTATCG ATCCACTGG	AAGGAGGAGGACGAGGATC	6377	6552	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	185	ACTAAGACTGTATCG ATCCACTGG	CCGCGCCTGAAGGAGGAG	6377	6561	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	191	ACTAAGACTGTATCG ATCCACTGG	TAGAGGCCCGCGCC TGAAG	6377	6567	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	106	CTAAGACTGTATCGA TCCACTGG	TCACGGGATTAGGCCCAATG	6378	6483	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	146	CTAAGACTGTATCGA TCCACTGG	ATGATCTGCACGG GTTGGG	6378	6523	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	160	CTAAGACTGTATCGA TCCACTGG	GAGTCCTGGGGA GCATGATC	6378	6537	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	181	CTAAGACTGTATCGA TCCACTGG	CGCCTGAAGGAG GAGGAC	6378	6558	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	186	ACATGGGGACTAAG ACTGTATCG	TGAAGGAGGAGG ACGAGGAG	6535	6554	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	178	ACTAAGACTGTATCG ATCCACTGG	TGAAGGAGGAGG ACGAGGAG	6535	6554	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	180	GGACTAAGACTGTA TCGATCCACTG	TGAAGGAGGAGG ACGAGGAG	6535	6554	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	184	ATGGGGACTAAGAC TGTATCGATC	TGAAGGAGGAGG ACGAGGAG	6535	6554	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	182	GGGGACTAAGACTG TATCGATCC	TGAAGGAGGAGG ACGAGGAG	6535	6554	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	191	AAACATGGGGACTA AGACTGTATCG	GCCTGAAGGAGG AGGACGAG	6538	6557	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	183	GGACTAAGACTGTA TCGATCCACTG	GCCTGAAGGAGG AGGACGAG	6538	6557	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	184	ACATGGGGACTAAG ACTGTATCG	GCCTGAAGGAGG AGGACGAG	6538	6557	no
XRMV env	76.9	TGACCCCTGTTCTCTCTGACCCGCCA	68.1	60	181	ACTAAGACTGTATCG ATCCACTGG	GCCTGAAGGAGG AGGACGAG	6538	6557	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.9	TGACCCTGTTCTCTCTGACCCGCCA	68.1	60	186	TGGGGACTAAGACTGTATCGATCC	GCCTGAAGGAGGAGGACGAG	6538	6557	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	95	CTTCAGGCGCGGCCTCTATG	GTAGGCTCCTTCTACGAGTTTAG	6550	6644	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	111	CTTCAGGCGCGGCCTCTATG	TGAGGTTGAGGGCTAGGTAGG	6550	6660	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	116	CTTCAGGCGCGGCCTCTATG	ACTGGTGAGGTGAGGGCTAG	6550	6665	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	121	CTTCAGGCGCGGCCTCTATG	TCGGGACTGGTGAAGGTGAG	6550	6670	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	123	CTTCAGGCGCGGCCTCTATG	TGTCGGGACTGGTGAGGTG	6550	6672	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	128	CTTCAGGCGCGGCCTCTATG	GGTTTTGTCGGGACTGGTGAG	6550	6677	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	129	CTTCAGGCGCGGCCTCTATG	GGGTTTTGTCGGGACTGGTG	6550	6678	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	133	CTTCAGGCGCGGCCTCTATG	TCTTGGGTTTTGTCGGGACTG	6550	6682	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	142	CTTCAGGCGCGGCCTCTATG	AGCCAGCACTCTTGGGTTTTG	6550	6691	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	150	CTTCAGGCGCGGCCTCTATG	CTAGACACAGCCAGCACTCTTG	6550	6699	no
XRMV env	77.1	CGCCTTCTCAACAACCTGGGACGGG	68.5	64	154	CTTCAGGCGCGGCCTCTATG	GATACTAGACACAGCCAGCACTC	6550	6703	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	155	CTTCAGGCGCGGCCTCTATG	CGATACTAGACAC AGCCAGCAC	6550	6704	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	89	TCAGGCGCGGCCTCTATG	GCTCCTTCTACCAG GTTTAGCAG	6552	6640	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	93	TCAGGCGCGGCCTCTATG	GTAGGCTCCTTCT ACCAGGTTAG	6552	6644	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	100	TCAGGCGCGGCCTCTATG	GGGCTAGGTAGG CTCCTTCTAC	6552	6651	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	97	GGCGCGGCCTCTATGGTG	GGGCTAGGTAGG CTCCTTCTAC	6552	6651	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	101	TCAGGCGCGGCCTCTATG	AGGGCTAGGTAG GCTCCTTC	6552	6652	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	109	TCAGGCGCGGCCTCTATG	TGAGGTTAGGGC TAGGTAGG	6552	6660	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	114	TCAGGCGCGGCCTCTATG	ACTGGTGAGGTTG AGGGCTAG	6552	6665	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	138	TCAGGCGCGGCCTCTATG	CCAGCACTCTTGG GTTTGTGTC	6552	6669	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	119	TCAGGCGCGGCCTCTATG	TCGGGACTGGTGA GGTTGAG	6552	6670	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	121	TCAGGCGCGGCCTCTATG	TGTCGGGACTGGT GAGGTTG	6552	6672	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	126	TCAGGCGCGGCCTCTATG	GGTTTTGTCGGGA CTGGTGAG	6552	6677	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	127	TCAGGCGCGGCCTCT ATG	GGGTTTGTGCGG ACTGGTG	6552	6678	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	131	TCAGGCGCGGCCTCT ATG	TCTTGGGTTTGTG GGGACTG	6552	6682	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	140	TCAGGCGCGGCCTCT ATG	AGCCAGCACTCTT GGGTTTGTG	6552	6691	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	143	TCAGGCGCGGCCTCT ATG	CACAGCCAGCACT CTTGGG	6552	6694	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	148	TCAGGCGCGGCCTCT ATG	CTAGACACAGCCA GCACTCTTG	6552	6699	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	152	TCAGGCGCGGCCTCT ATG	GATACTAGACACA GCCAGCACTC	6552	6703	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	153	TCAGGCGCGGCCTCT ATG	CGATACTAGACAC AGCCAGCAC	6552	6704	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	159	TCAGGCGCGGCCTCT ATG	GGGTCCCAGTACT AGACACAG	6552	6710	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	161	TCAGGCGCGGCCTCT ATG	GGGGTCCCGATA CTAGACAC	6552	6712	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	182	TTCAGGCGCGGCCTC TATG	CGGCCACCCCTTC GTAGTAG	6552	6732	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	187	TCAGGCGCGGCCTCT ATG	CTAGGACGGCCAC CCCTTC	6552	6738	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	180	CAGGCGCGGCCTCT ATGG	CGGCCACCCCTTC GTAGTAG	6553	6732	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	86	GGCGCGGCCTCTAT GGTG	GCTCCTTCTACCAG GTTTAGCAG	6555	6640	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	90	GGCGCGGCCTCTAT GGTG	GTAGGCTCCTTCT ACCAGGTTTAG	6555	6644	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	106	GGCGCGGCCTCTAT GGTG	TGAGGTTGAGGGC TAGGTAGG	6555	6660	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	111	GGCGCGGCCTCTAT GGTG	ACTGGTGAGGTG AGGGCTAG	6555	6665	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	135	GGCGCGGCCTCTAT GGTG	CCAGCACTCTTG GTTTTGTC	6555	6669	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	116	GGCGCGGCCTCTAT GGTG	TCGGGACTGGTGA GGTTGAG	6555	6670	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	118	GGCGCGGCCTCTAT GGTG	TGTCGGGACTGGT GAGGTTG	6555	6672	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	123	GGCGCGGCCTCTAT GGTG	GGTTTTGTCGGGA CTGGTGAG	6555	6677	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	124	GGCGCGGCCTCTAT GGTG	GGGTTTTGTCGGG ACTGGTG	6555	6678	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	128	GGCGCGGCCTCTAT GGTG	TCTTGGGTTTTGTC GGGACTG	6555	6682	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	137	GGCGCGGCCTCTAT GGTG	AGCCAGCACTCTT GGGTTTTG	6555	6691	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	140	GGCGCGGCCTCTAT GGTG	CACAGCCAGCACT CTTGGG	6555	6694	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	145	GGCGGGCCTCTAT GGTG	CTAGACACAGCCA GCACTCTTG	6555	6699	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	149	GGCGGGCCTCTAT GGTG	GATACTAGACACA GCCAGCACTC	6555	6703	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	150	GGCGGGCCTCTAT GGTG	CGATACTAGACAC AGCCAGCAC	6555	6704	no
XRMV env	77.1	CGCCTTCTCAACA ACCTGGGACGGG	68.5	64	178	GGCGGGCCTCTAT GGTG	CGGCCACCCCTTC GTAGTAG	6555	6732	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	159	GACGGGAGACAGGC TGCTAAAC	GCAGAGGTATGGT TGGAGTAAGTAC	6605	6763	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	159	GACGGGAGACAGGC TGCTAAAC	GCAGAGGTATGGT TGGAGTAAGTAC	6605	6763	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	127	ACGGGAGACAGGCT GCTAAAC	CGGCCACCCCTTC GTAGTAG	6606	6732	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	93	CGGGAGACAGGCTG CTAAAC	CTAGACACAGCCA GCACTCTTG	6607	6699	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	93	CGGGAGACAGGCTG CTAAAC	CTAGACACAGCCA GCACTCTTG	6607	6699	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	97	CGGGAGACAGGCTG CTAAAC	GATACTAGACACA GCCAGCACTC	6607	6703	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	97	CGGGAGACAGGCTG CTAAAC	GATACTAGACACA GCCAGCACTC	6607	6703	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	127	ACGGGAGACAGGCT GCTAAAC	CGGCCACCCCTTC GTAGTAG	6607	6732	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	157	CGGGAGACAGGCTG CTAAAC	GCAGAGGTATGGT TGGAGTAAGTAC	6607	6763	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	157	CGGGAGACAGGCTG CTAAAC	GCAGAGGTATGGT TGGAGTAAGTAC	6607	6763	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	158	CGGGAGACAGGCTG CTAAAC	GGCAGAGGTATG GTTGGAGTAAG	6607	6764	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	158	CGGGAGACAGGCTG CTAAAC	GGCAGAGGTATG GTTGGAGTAAG	6607	6764	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	159	CGGGAGACAGGCTG CTAAAC	GGCAGAGGTAT GGTTGGAG	6607	6765	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	159	CGGGAGACAGGCTG CTAAAC	GGCAGAGGTAT GGTTGGAG	6607	6765	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	161	CGGGAGACAGGCTG CTAAAC	CGGGCAGAGGT ATGGTTG	6607	6767	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	161	CGGGAGACAGGCTG CTAAAC	CGGGCAGAGGT ATGGTTG	6607	6767	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	163	CGGGAGACAGGCTG CTAAAC	GCCGGGCAGAG GTATGG	6607	6769	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	163	CGGGAGACAGGCTG CTAAAC	GCCGGGCAGAG GTATGG	6607	6769	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	185	CGGGAGACAGGCTG CTAAAC	TTGGGAGGTCAG GAGCAG	6607	6791	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	185	CGGGAGACAGGCTG CTAAAC	TTGGGAGGTCAG GAGCAG	6607	6791	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	192	CGGAGACAGGCTG CTAAAC	GCTTGTGTTGGGA GGTCACG	6607	6798	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	192	CGGAGACAGGCTG CTAAAC	GCTTGTGTTGGGA GGTCACG	6607	6798	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	196	CGGAGACAGGCTG CTAAAC	GTCAAGCTTGTGTT GGGAGGTC	6607	6802	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	196	CGGAGACAGGCTG CTAAAC	GTCAAGCTTGTGTT GGGAGGTC	6607	6802	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	200	CGGAGACAGGCTG CTAAAC	CAGGGTCAGCTTG TGTTGGG	6607	6806	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	200	CGGAGACAGGCTG CTAAAC	CAGGGTCAGCTTG TGTTGGG	6607	6806	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	103	GGGAGACAGGCTGC TAAACC	GGGTCCCATACT AGACACAG	6608	6710	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	103	GGGAGACAGGCTGC TAAACC	GGGTCCCATACT AGACACAG	6608	6710	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	178	GGGAGACAGGCTGC TAAACC	GGTCACGGAGCA GTTAGCC	6608	6785	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	178	GGGAGACAGGCTGC TAAACC	GGTCACGGAGCA GTTAGCC	6608	6785	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	181	GGGAGACAGGCTGC TAAACC	GGAGGTCACGGA GCAGTTAG	6608	6788	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	181	GGGAGACAGGCTGC TAAACC	GGAGGTCACGGA GCAGTTAG	6608	6788	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	152	GACAGGCTGCTAAA CCTGGTAG	GCAGAGGTATGGT TGGAGTAAGTAC	6612	6763	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	152	GACAGGCTGCTAAA CCTGGTAG	GCAGAGGTATGGT TGGAGTAAGTAC	6612	6763	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	153	GACAGGCTGCTAAA CCTGGTAG	GGCAGAGGTATG GTTGGAGTAAG	6612	6764	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	153	GACAGGCTGCTAAA CCTGGTAG	GGCAGAGGTATG GTTGGAGTAAG	6612	6764	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	195	GACAGGCTGCTAAA CCTGGTAG	CAGGGTCAGCTTG TGTTGGG	6612	6806	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	195	GACAGGCTGCTAAA CCTGGTAG	CAGGGTCAGCTTG TGTTGGG	6612	6806	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	98	ACAGGCTGCTAAAC CTGGTAG	GGGTCCCATACT AGACACAG	6613	6710	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	98	ACAGGCTGCTAAAC CTGGTAG	GGGTCCCATACT AGACACAG	6613	6710	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	153	ACAGGCTGCTAAAC CTGGTAG	GGGCAGAGGTAT GGTTGGAG	6613	6765	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	153	ACAGGCTGCTAAAC CTGGTAG	GGGCAGAGGTAT GGTTGGAG	6613	6765	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	176	ACAGGCTGCTAAAC CTGGTAG	GGAGGTCACGGA GCAGTTAG	6613	6788	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	176	ACAGGCTGCTAAAC CTGGTAG	GGAGGTCACGGA GCAGTTAG	6613	6788	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	85	AGGCTGCTAAACCT GGTAGAAG	CTAGACACAGCCA GCACTCTTG	6615	6699	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	85	AGGCTGCTAAACCT GGTAGAAG	CTAGACACAGCCA GCACTCTTG	6615	6699	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	89	AGGCTGCTAAACCT GGTAGAAG	GATACTAGACACA GCCAGCACTC	6615	6703	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	89	AGGCTGCTAAACCT GGTAGAAG	GATACTAGACACA GCCAGCACTC	6615	6703	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	96	AGGCTGCTAAACCT GGTAGAAG	GGGTCCCATACT AGACACAG	6615	6710	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	96	AGGCTGCTAAACCT GGTAGAAG	GGGTCCCATACT AGACACAG	6615	6710	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	150	CAGGCTGCTAAACCT GGTAGAAG	GCAGAGGTATGGT TGGAGTAAGTAC	6615	6763	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	150	CAGGCTGCTAAACCT GGTAGAAG	GCAGAGGTATGGT TGGAGTAAGTAC	6615	6763	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	150	AGGCTGCTAAACCT GGTAGAAG	GGCAGAGGTATG GTTGGAGTAAG	6615	6764	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	150	AGGCTGCTAAACCT GGTAGAAG	GGCAGAGGTATG GTTGGAGTAAG	6615	6764	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	151	AGGCTGCTAAACCT GGTAGAAG	GGGCAGAGGTAT GGTTGGAG	6615	6765	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	151	AGGCTGCTAAACCT GGTAGAAG	GGGCAGAGGTAT GGTTGGAG	6615	6765	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	153	AGGCTGCTAAACCT GGTAGAAG	CGGGCAGAGGT ATGTTG	6615	6767	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	153	AGGCTGCTAAACCT GGTAGAAG	CGGGCAGAGGT ATGTTG	6615	6767	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	155	AGGCTGCTAAACCT GGTAGAAG	GCCGGGCAGAG GTATGG	6615	6769	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	155	AGGCTGCTAAACCT GGTAGAAG	GCCGGGCAGAG GTATGG	6615	6769	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	171	AGGCTGCTAAACCT GGTAGAAG	GGTCACGGAGCA GTTAGCC	6615	6785	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	171	AGGCTGCTAAACCT GGTAGAAG	GGTCACGGAGCA GTTAGCC	6615	6785	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	174	AGGCTGCTAAACCT GGTAGAAG	GGAGGTCACGGA GCAGTTAG	6615	6788	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	174	AGGCTGCTAAACCT GGTAGAAG	GGAGGTCACGGA GCAGTTAG	6615	6788	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	184	AGGCTGCTAAACCT GGTAGAAG	GCTTGTGTTGGGA GGTCACG	6615	6798	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	184	AGGCTGCTAAACCT GGTAGAAG	GCTTGTGTTGGGA GGTCACG	6615	6798	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	188	AGGCTGCTAAACCT GGTAGAAG	GTCAGCTTGTGTT GGGAGGTC	6615	6802	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	188	AGGCTGCTAAACCT GGTAGAAG	GTCAGCTTGTGTT GGGAGGTC	6615	6802	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	192	AGGCTGCTAAACCT GGTAGAAG	CAGGGTCAGCTTG TGTTGGG	6615	6806	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	192	AGGCTGCTAAACCT GGTAGAAG	CAGGGTCAGCTTG TGTTGGG	6615	6806	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	195	AGGCTGCTAAACCT GGTAGAAG	GGACAGGGTCAG CTTGTGTG	6615	6809	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	195	AGGCTGCTAAACCT GGTAGAAG	GGACAGGGTCAG CTTGTGTG	6615	6809	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	196	AGGCTGCTAAACCT GGTAGAAG	CGGACAGGGTCA GCTTGTG	6615	6810	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	196	AGGCTGCTAAACCT GGTAGAAG	CGGACAGGGTCA GCTTGTG	6615	6810	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	148	GGCTGCTAAACCTG GTAGAAGG	GCAGAGGTATGGT TGGAGTAAGTAC	6616	6763	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	148	GGCTGCTAAACCTG GTAGAAGG	GCAGAGGTATGGT TGGAGTAAGTAC	6616	6763	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	83	GCTGCTAAACCTGGT AGAAGGAG	CTAGACACAGCCA GCACTCTTG	6617	6699	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	83	GCTGCTAAACCTGGT AGAAGGAG	CTAGACACAGCCA GCACTCTTG	6617	6699	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	87	GCTGCTAAACCTGGT AGAAGGAG	GATACTAGACACA GCCAGCACTC	6617	6703	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	87	GCTGCTAAACCTGGT AGAAGGAG	GATACTAGACACA GCCAGCACTC	6617	6703	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	148	GCTGCTAAACCTGGT AGAAGGAG	GGCAGAGGTATG GTTGGAGTAAG	6617	6764	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	148	GCTGCTAAACCTGGT AGAAGGAG	GGCAGAGGTATG GTTGGAGTAAG	6617	6764	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	190	GCTGCTAAACCTGGT AGAAGGAG	CAGGGTCAGCTTG TGTTGGG	6617	6806	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	190	GCTGCTAAACCTGGT AGAAGGAG	CAGGGTCAGCTTG TGTTGGG	6617	6806	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	193	GCTGCTAAACCTGGT AGAAGGAG	GGACAGGGTCAG CTTGTGTG	6617	6809	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	193	GCTGCTAAACCTGGT AGAAGGAG	GGACAGGGTCAG CTTGTGTG	6617	6809	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	79	CTAAACCTGGTAGA AGGAGCCTAC	CTAGACACAGCCA GCACTCTTG	6621	6699	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	79	CTAAACCTGGTAGA AGGAGCCTAC	CTAGACACAGCCA GCACTCTTG	6621	6699	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	83	CTAAACCTGGTAGA AGGAGCCTAC	GATACTAGACACA GCCAGCACTC	6621	6703	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	83	CTAAACCTGGTAGA AGGAGCCTAC	GATACTAGACACA GCCAGCACTC	6621	6703	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	143	CTAAACCTGGTAGA AGGAGCCTAC	GCAGAGGTATGGT TGGAGTAAGTAC	6621	6763	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	143	CTAAACCTGGTAGA AGGAGCCTAC	GCAGAGGTATGGT TGGAGTAAGTAC	6621	6763	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	144	CTAAACCTGGTAGA AGGAGCCTAC	GGCAGAGGTATG GTTGGAGTAAG	6621	6764	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	144	CTAAACCTGGTAGA AGGAGCCTAC	GGCAGAGGTATG GTTGGAGTAAG	6621	6764	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	165	CTAAACCTGGTAGA AGGAGCCTAC	GGTCACGGAGCA GTTAGGCC	6621	6785	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	165	CTAAACCTGGTAGA AGGAGCCTAC	GGTCACGGAGCA GTTAGGCC	6621	6785	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	167	TAAACCTGGTAGAA GGAGCCTAC	GGAGGTCACGGA GCAGTTAG	6622	6788	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	167	TAAACCTGGTAGAA GGAGCCTAC	GGAGGTCACGGA GCAGTTAG	6622	6788	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	163	AAACCTGGTAGAAG GAGCCTAC	GGTCACGGAGCA GTTAGGCC	6623	6785	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	163	AAACCTGGTAGAAG GAGCCTAC	GGTCACGGAGCA GTTAGGCC	6623	6785	no
XRMV env	76.8	AGCCCTCAACCTC ACCAGTCCCGA	68	62.5	165	AACCTGGTAGAAGG AGCCTAC	GGAGGTCACGGA GCAGTTAG	6624	6788	no
XRMV env	76.5	TAGCCCTCAACCTC ACCAGTCCCGA	67.7	60	165	AACCTGGTAGAAGG AGCCTAC	GGAGGTCACGGA GCAGTTAG	6624	6788	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	88	GCATAGGAGCAGTT CCCAAAAC	GCGGGAGAGGCC AAATAGTAG	6832	6919	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	88	GCATAGGAGCAGTT CCCAAAAC	GCGGGAGAGGCC AAATAGTAG	6832	6919	no
XRMV env	75.9	ACCACCCAGAAGA	67.2	65.2	113	GCATAGGAGCAGTT	GGTGCTGCAAGCC	6832	6944	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
		CGAGCGACGG				CCCCAAAC	CAATG			
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	113	GCATAGGAGCAGTT CCCCAAAC	GGTGTGCAAGCC CAATG	6832	6944	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	141	GCATAGGAGCAGTT CCCCAAAC	CAGTAGTAGATAG ACAGGGAGTGAG	6832	6972	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	141	GCATAGGAGCAGTT CCCCAAAC	CAGTAGTAGATAG ACAGGGAGTGAG	6832	6972	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	160	GCATAGGAGCAGTT CCCCAAAC	TCAGTGGTTAAGT TAAGCACAGTAG	6832	6991	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	160	GCATAGGAGCAGTT CCCCAAAC	TCAGTGGTTAAGT TAAGCACAGTAG	6832	6991	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	173	GCATAGGAGCAGTT CCCCAAAC	CAGGACACAGTAA TCAGTGGTTAAG	6832	7004	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	173	GCATAGGAGCAGTT CCCCAAAC	CAGGACACAGTAA TCAGTGGTTAAG	6832	7004	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	178	GCATAGGAGCAGTT CCCCAAAC	TCAACCCAGGACAC AGTAATCAGTG	6832	7009	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	178	GCATAGGAGCAGTT CCCCAAAC	TCAACCCAGGACAC AGTAATCAGTG	6832	7009	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	180	GCATAGGAGCAGTT CCCCAAAC	GTTCAACCCAGGAC ACAGTAATCAG	6832	7011	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	180	GCATAGGAGCAGTT CCCCAAAC	GTTCAACCCAGGAC ACAGTAATCAG	6832	7011	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	174	GAGCAGTTCCTCAA ACCCATC	GTTCAACCCAGGAC ACAGTAATCAG	6832	7011	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	182	GCATAGGAGCAGTT CCCCAAAC	GAGTTCAACCAGG ACACAGTAATC	6832	7013	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	182	GCATAGGAGCAGTT CCCCAAAC	GAGTTCAACCAGG ACACAGTAATC	6832	7013	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	185	GCATAGGAGCAGTT CCCCAAAC	CCAGAGTTCAACC AGGACACAG	6832	7016	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	185	GCATAGGAGCAGTT CCCCAAAC	CCAGAGTTCAACC AGGACACAG	6832	7016	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	186	GCATAGGAGCAGTT CCCCAAAC	GCCAGAGTTCAAC CAGGACAC	6832	7017	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	186	GCATAGGAGCAGTT CCCCAAAC	GCCAGAGTTCAAC CAGGACAC	6832	7017	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	187	GCATAGGAGCAGTT CCCCAAAC	GGCCAGAGTTCAA CCAGGAC	6832	7018	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	187	GCATAGGAGCAGTT CCCCAAAC	GGCCAGAGTTCAA CCAGGAC	6832	7018	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	147	CATAGGAGCAGTTC CCAAACCC	TTAAGCACAGTAG TAGATAGACAGG	6833	6979	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	147	CATAGGAGCAGTTC CCAAACCC	TTAAGCACAGTAG TAGATAGACAGG	6833	6979	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	171	ATAGGAGCAGTTCC CAAAACCC	CAGGACACAGTAA TCAGTGGTTAAG	6834	7004	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	171	ATAGGAGCAGTTCC CAAAACCC	CAGGACACAGTAA TCAGTGGTTAAG	6834	7004	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	176	ATAGGAGCAGTTCC CAAAACCC	TCAACCAGGACAC AGTAATCAGTG	6834	7009	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	176	ATAGGAGCAGTTCC CAAAACCC	TCAACCAGGACAC AGTAATCAGTG	6834	7009	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	178	ATAGGAGCAGTTCC CAAAACCC	GTTCAACCAGGAC ACAGTAATCAG	6834	7011	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	178	ATAGGAGCAGTTCC CAAAACCC	GTTCAACCAGGAC ACAGTAATCAG	6834	7011	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	183	ATAGGAGCAGTTCC CAAAACCC	CCAGAGTTCAACC AGGACACAG	6834	7016	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	183	ATAGGAGCAGTTCC CAAAACCC	CCAGAGTTCAACC AGGACACAG	6834	7016	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	184	ATAGGAGCAGTTCC CAAAACCC	GCCAGAGTTCAAC CAGGACAC	6834	7017	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	184	ATAGGAGCAGTTCC CAAAACCC	GCCAGAGTTCAAC CAGGACAC	6834	7017	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	185	ATAGGAGCAGTTCC CAAAACCC	GGCCAGAGTTCAA CCAGGAC	6834	7018	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	185	ATAGGAGCAGTTCC CAAAACCC	GGCCAGAGTTCAA CCAGGAC	6834	7018	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	157	TAGGAGCAGTTCCC AAAACCC	TCAGTGGTTAAGT TAAGCACAGTAG	6835	6991	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	157	TAGGAGCAGTTCCC AAAACCC	TCAGTGGTTAAGT TAAGCACAGTAG	6836	6991	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	75.9	ACCACCCAGAAAGACGAGCGACGG	67.2	65.2	169	AGGAGCAGTCCCAAAACCC	CAGGACACAGTAA TCAGTGGTTAAG	6836	7004	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	169	AGGAGCAGTCCCAAAACCC	CAGGACACAGTAA TCAGTGGTTAAG	6836	7004	no
XRMV env	75.9	ACCACCCAGAAAGACGAGCGACGG	67.2	65.2	176	AGGAGCAGTCCCAAAACCC	GTTCAACCCAGGAC ACAGTAATCAG	6836	7011	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	176	AGGAGCAGTCCCAAAACCC	GTTCAACCCAGGAC ACAGTAATCAG	6836	7011	no
XRMV env	75.9	ACCACCCAGAAAGACGAGCGACGG	67.2	65.2	183	AGGAGCAGTCCCAAAACCC	GGCCAGAGTTCAA CCAGGAC	6836	7018	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	183	AGGAGCAGTCCCAAAACCC	GGCCAGAGTTCAA CCAGGAC	6836	7018	no
XRMV env	75.9	ACCACCCAGAAAGACGAGCGACGG	67.2	65.2	107	GAGCAGTCCCAAA ACCCATC	GGTGTGCAAGCC CAAATG	6838	6944	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	107	GAGCAGTCCCAAA ACCCATC	GGTGTGCAAGCC CAAATG	6838	6944	no
XRMV env	75.9	ACCACCCAGAAAGACGAGCGACGG	67.2	65.2	135	GAGCAGTCCCAAA ACCCATC	CAGTAGTAGATAG ACAGGGAGTGAG	6838	6972	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	135	GAGCAGTCCCAAA ACCCATC	CAGTAGTAGATAG ACAGGGAGTGAG	6838	6972	no
XRMV env	75.9	ACCACCCAGAAAGACGAGCGACGG	67.2	65.2	142	GAGCAGTCCCAAA ACCCATC	TTAAGCACAGTAG TAGATAGACAGG	6838	6979	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	142	GAGCAGTCCCAAA ACCCATC	TTAAGCACAGTAG TAGATAGACAGG	6838	6979	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	75.9	ACCACCCAGAAAGACGACGG	67.2	65.2	151	GAGCAGTTCCTCAAAACCCATC	GTGGTTAAGTTAAGCACAGTAGTAG	6838	6988	no
XRMV env	75.7	ATACCACCCAGAAAGACGACGG	67.1	60	151	GAGCAGTTCCTCAAAACCCATC	GTGGTTAAGTTAAGCACAGTAGTAG	6838	6988	no
XRMV env	75.9	ACCACCCAGAAAGACGACGG	67.2	65.2	154	GAGCAGTTCCTCAAAACCCATC	TCAGTGGTTAAGTTAAGCACAGTAG	6838	6991	no
XRMV env	75.7	ATACCACCCAGAAAGACGACGG	67.1	60	154	GAGCAGTTCCTCAAAACCCATC	TCAGTGGTTAAGTTAAGCACAGTAG	6838	6991	no
XRMV env	75.9	ACCACCCAGAAAGACGACGG	67.2	65.2	161	GAGCAGTTCCTCAAAACCCATC	ACAGTAATCAGTGTTAAGTTAAGC	6838	6998	no
XRMV env	75.7	ATACCACCCAGAAAGACGACGG	67.1	60	161	GAGCAGTTCCTCAAAACCCATC	ACAGTAATCAGTGTTAAGTTAAGC	6838	6998	no
XRMV env	75.9	ACCACCCAGAAAGACGACGG	67.2	65.2	167	GAGCAGTTCCTCAAAACCCATC	CAGGACACAGTAATCAGTGTTAAG	6838	7004	no
XRMV env	75.7	ATACCACCCAGAAAGACGACGG	67.1	60	167	GAGCAGTTCCTCAAAACCCATC	CAGGACACAGTAATCAGTGTTAAG	6838	7004	no
XRMV env	75.9	ACCACCCAGAAAGACGACGG	67.2	65.2	173	GGAGCAGTTCCTCAAAACCCATC	TCAACCCAGGACACAGTAATCAGTG	6838	7009	no
XRMV env	75.7	ATACCACCCAGAAAGACGACGG	67.1	60	173	GGAGCAGTTCCTCAAAACCCATC	TCAACCCAGGACACAGTAATCAGTG	6838	7009	no
XRMV env	75.9	ACCACCCAGAAAGACGACGG	67.2	65.2	174	GAGCAGTTCCTCAAAACCCATC	GTTCAACCCAGGACACAGTAATCAG	6838	7011	no
XRMV env	75.9	ACCACCCAGAAAGACGACGG	67.2	65.2	176	GAGCAGTTCCTCAAAACCCATC	GAGTTCAACCCAGGACACAGTAATC	6838	7013	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	176	GAGCAGTTCCCAA ACCCATC	GAGTTCAACCAGG ACACAGTAATC	6838	7013	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	180	GGAGCAGTTCCCAA AACCACATC	CCAGAGTTCAACC AGGACACAG	6838	7016	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	180	GGAGCAGTTCCCAA AACCACATC	CCAGAGTTCAACC AGGACACAG	6838	7016	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	181	GGAGCAGTTCCCAA AACCACATC	GCCAGAGTTCAAC CAGGACAC	6838	7017	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	181	GAGCAGTTCCCAA ACCCATC	GGCCAGAGTTCAA CCAGGAC	6838	7018	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	181	GAGCAGTTCCCAA ACCCATC	GGCCAGAGTTCAA CCAGGAC	6838	7018	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	171	AGCAGTTCCCAAAC CCATCAG	TCAACCCAGGACAC AGTAATCAGTG	6839	7009	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	80	GCAGTTCCCAAAC CATCAG	GCGGGAGAGGCC AAATAGTAG	6840	6919	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	80	GCAGTTCCCAAAC CATCAG	GCGGGAGAGGCC AAATAGTAG	6840	6919	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	81	GCAGTTCCCAAAC CATCAG	GCGGGAGAGGCC CAAATAG	6840	6920	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	81	GCAGTTCCCAAAC CATCAG	GCGGGAGAGGCC CAAATAG	6840	6920	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	105	GCAGTTCCCAAAC CATCAG	GGTGTGCAAGCC CAAATG	6840	6944	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	105	GCAGTTCCTCCAAAACCC CATCAG	GGTGCTGCAAGCC CAATG	6840	6944	no
XRMV env	75.9	ACCACCCAGAAAGA CGAGCGACGG	67.2	65.2	133	GCAGTTCCTCCAAAACCC CATCAG	CAGTAGTAGATAG ACAGGGAGTGAG	6840	6972	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	133	GCAGTTCCTCCAAAACCC CATCAG	CAGTAGTAGATAG ACAGGGAGTGAG	6840	6972	no
XRMV env	75.9	ACCACCCAGAAAGA CGAGCGACGG	67.2	65.2	140	GCAGTTCCTCCAAAACCC CATCAG	TTAAGCACAGTAG TAGATAGACAGG	6840	6979	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	140	GCAGTTCCTCCAAAACCC CATCAG	TTAAGCACAGTAG TAGATAGACAGG	6840	6979	no
XRMV env	75.9	ACCACCCAGAAAGA CGAGCGACGG	67.2	65.2	149	GCAGTTCCTCCAAAACCC CATCAG	GTGGTTAAGTTAA GCACAGTAGTAG	6840	6988	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	149	GCAGTTCCTCCAAAACCC CATCAG	GTGGTTAAGTTAA GCACAGTAGTAG	6840	6988	no
XRMV env	75.9	ACCACCCAGAAAGA CGAGCGACGG	67.2	65.2	152	GCAGTTCCTCCAAAACCC CATCAG	TCAGTGGTTAAGT TAAGCACAGTAG	6840	6991	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	152	GCAGTTCCTCCAAAACCC CATCAG	TCAGTGGTTAAGT TAAGCACAGTAG	6840	6991	no
XRMV env	75.9	ACCACCCAGAAAGA CGAGCGACGG	67.2	65.2	165	GCAGTTCCTCCAAAACCC CATCAG	CAGGACACAGTAA TCAGTGGTTAAG	6840	7004	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	165	GCAGTTCCTCCAAAACCC CATCAG	CAGGACACAGTAA TCAGTGGTTAAG	6840	7004	no
XRMV env	75.7	ATACCAACCCAGAA GACGAGCGACGG	67.1	60	171	AGCAGTTCCTCCAAAACCC CCATCAG	TCAACCCAGGACAC AGTAATCAGTG	6840	7009	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	172	GCAGTTCCTCCAAAACC CATCAG	GTTCAACCAGGAC ACAGTAATCAG	6840	7011	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	172	GCAGTTCCTCCAAAACC CATCAG	GTTCAACCAGGAC ACAGTAATCAG	6840	7011	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	174	GCAGTTCCTCCAAAACC CATCAG	GAGTTCAACCAGG ACACAGTAATC	6840	7013	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	174	GCAGTTCCTCCAAAACC CATCAG	GAGTTCAACCAGG ACACAGTAATC	6840	7013	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	178	AGCAGTTCCTCCAAAACC CCATCAG	CCAGAGTTCAACC AGGACACAG	6840	7016	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	178	AGCAGTTCCTCCAAAACC CCATCAG	CCAGAGTTCAACC AGGACACAG	6840	7016	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	179	AGCAGTTCCTCCAAAACC CCATCAG	GCCAGAGTTCAAC CAGGACAC	6840	7017	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	181	GGAGCAGTTCCTCCAA AACCCTATC	GCCAGAGTTCAAC CAGGACAC	6840	7017	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	179	AGCAGTTCCTCCAAAACC CCATCAG	GCCAGAGTTCAAC CAGGACAC	6840	7017	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	179	GCAGTTCCTCCAAAACC CATCAG	GGCCAGAGTTCAA CCAGGAC	6840	7018	no
XRMV env	75.7	ATACCAACCAGAA GACGAGCGACGG	67.1	60	179	GCAGTTCCTCCAAAACC CATCAG	GGCCAGAGTTCAA CCAGGAC	6840	7018	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	158	CAGTTCCTCCAAAACC ATCAGG	ACAGTAATCAGTG GTTAAGTTAAGC	6841	6998	no

(continued)

Sequence	Rating	Probe Sequence	T _m	GC%	Product Length	Sense Primer	Anti-sense Primer	5' XMRV coordinate	3' XMRV coordinate	Detects MLV related and XMRV
XRMV env	75.7	ATACCACCCAGAA GACGAGCGACGG	67.1	60	158	CAGTTCCCAAAACCC ATCAGG	ACAGTAATCAGTG GTTAAGTTAAGC	6841	6998	no
XRMV env	75.9	ACCACCCAGAAGA CGAGCGACGG	67.2	65.2	176	CAGTTCCCAAAACCC ATCAGG	CCAGAGTTCAACC AGGACACAG	6841	7016	no
XRMV env	75.7	ATACCACCCAGAA GACGAGCGACGG	67.1	60	176	CAGTTCCCAAAACCC ATCAGG	CCAGAGTTCAACC AGGACACAG	6841	7016	no

[0116] RT-PCR assay is performed to quantify RNA from biological samples, such as whole blood and plasma, without the need for RNA extraction. The assay employs a two-step amplification process with the initial step consisting of the distribution of 2 μ L of experimental sample directly into a cDNA reaction mix. Following completion of the reverse transcriptase (RT) cDNA synthesis, a 2 μ L aliquot is removed, transferred into a qPCR reaction mix and a qPCR protocol is performed.

[0117] In a separate aspect, RNA is isolated from cell culture supernatants, whole blood or plasma using QIAamp®viral RNA Mini kit. The RNA is then used for reverse transcription. After reverse transcription, cDNA is subjected to qPCR assay.

SEQUENCE LISTING

[0118]

<110> Tocagen Inc.
Jolly, Douglas
Lin, Amy
Perez, Omar

<120> RETROVIRUS DETECTION

<130> 00014-010WO1

<140> Not yet known

<141> 2011-07-16

<150> US 61/365,297

<151> 2010-07-16

<150> 61/386,941

<151> 2010-09-27

<150> 61/391,360

<151> 2010-10-08

<160> 25

<170> PatentIn version 3.5

<210> 1

<211> 21

<212> DNA

<213> Artificial Sequence

<220>

<223> Oligonucleotide primer

<400> 1

actaccacctc tgagtctaac c 21

<210> 2

<211> 22

<212> DNA

<213> Artificial Sequence

<220>

<223> Oligonucleotide primer

<400> 2

ggccatccta cattgaaagt tg 22

<210> 3
 <211> 24
 <212> DNA
 <213> Artificial Sequence
 5
 <220>
 <223> Oligonucleotide Probe
 <400> 3
 10 cgattgcat ccaaccagtc tgtg 24
 <210> 4
 <211> 20
 <212> DNA
 15 <213> Artificial Sequence
 <220>
 <223> Oligonucleotide primer
 20 <400> 4
 ttgattcct cagtgggctc 20
 <210> 5
 <211> 22
 25 <212> DNA
 <213> Artificial Sequence
 <220>
 <223> Oligonucleotide primer
 30 <400> 5
 cgatacagtc ttagtcccca tg 22
 <210> 6
 35 <211> 26
 <212> DNA
 <213> Artificial Sequence
 <220>
 40 <223> Oligonucleotide probe
 <400> 6
 cccttttacc cgcgtcagtg aattct 26
 45 <210> 7
 <211> 21
 <212> DNA
 <213> Artificial Sequence
 50 <220>
 <223> Oligonucleotide primer
 <400> 7
 55 aacaagcggg tggaagacat c 21
 <210> 8
 <211> 21
 <212> DNA

<213> Artificial Sequence
 <220>
 <223> Oligonucleotide primer
 5
 <400> 8
 caaaggcgaa gagaggctga c 21
 <210> 9
 10 <211> 25
 <212> DNA
 <213> Artificial Sequence
 <220>
 15 <223> Oligonucleotide probe
 <400> 9
 cccaccgtgc ccaacccta caacc 25
 20 <210> 10
 <211> 21
 <212> DNA
 <213> Artificial Sequence
 25 <220>
 <223> Oligonucleotide primer
 <400> 10
 30 caaggggcta ctggaggaaa g 21
 <210> 11
 <211> 22
 <212> DNA
 <213> Artificial Sequence
 35 <220>
 <223> Oligonucleotide primer
 <400> 11
 40 ctttctcca tgtaccagac tg 22
 <210> 12
 <211> 23
 <212> DNA
 45 <213> Artificial Sequence
 <220>
 <223> Oligonucleotide probe
 50 <400> 12
 tatcgctgga ccacggatcg caa 23
 <210> 13
 <211> 22
 55 <212> DNA
 <213> Artificial Sequence
 <220>

<223> Oligonucleotide primer

 <400> 13
 accctcaacc gccctacaa gt 22
 5

 <210> 14
 <211> 20
 <212> DNA
 <213> Artificial Sequence
 10

 <220>
 <223> Oligonucleotide primer

 <400> 14
 15 gttaagcgcc tgataggctc 20

 <210> 15
 <211> 26
 <212> DNA
 20 <213> Artificial Sequence

 <220>
 <223> Oligonucleotide probe

 <400> 15
 25 ccccaaatga aagacccccg ctgacg 26

 <210> 16
 <211> 19
 30 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> Oligonucleotide primer
 35

 <400> 16
 agcccacaac ccctcactc 19

 <210> 17
 40 <211> 18
 <212> DNA
 <213> Artificial Sequence

 <220>
 45 <223> Oligonucleotide primer

 <400> 17
 tctcccgatc ccggacga 18

 <210> 18
 50 <211> 26
 <212> DNA
 <213> Artificial Sequence

 <220>
 55 <223> Oligonucleotide probe

 <400> 18

ccccaaatga aagacccccg ctgacg 26

5 <210> 19
 <211> 23
 <212> DNA
 <213> Artificial Sequence

10 <220>
 <223> Oligonucleotide primer

<400> 19
 atcatcatgt acggcatccc tag 23

15 <210> 20
 <211> 24
 <212> DNA
 <213> Artificial Sequence

20 <220>
 <223> Oligonucleotide primer

<400> 20
 tgaactgctt catcagcttc ttac 24

25 <210> 21
 <211> 25
 <212> DNA
 <213> Artificial Sequence

30 <220>
 <223> Oligonucleotide probe

<400> 21
 tcatcgtcaa caaccaccac ctcgt 25

35 <210> 22
 <211> 11892
 <212> DNA
 <213> Artificial Sequence

40 <220>
 <223> RCR Vector - pAC3-yCD2

45 <400> 22

50

55

EP 2 913 402 B1

	tagttattaa tagtaatcaa ttacggggtc attagttcat agcccatata tggagttccg	60
	cgttacataa cttacggtaa atggcccgcc tggctgaccg cccaacgacc cccgccatt	120
5	gacgtcaata atgacgtatg ttcccatagt aacgccaata gggactttcc attgacgtca	180
	atgggtggag tatttacggt aaactgccc cttggcagta catcaagtgt atcatatgcc	240
	aagtacgccc cctattgacg tcaatgacgg taaatggccc gcctggcatt atgccagta	300
10	catgacctta tgggactttc ctacttggca gtacatctac gtattagtca tcgctattac	360
	catggtgatg cgggttttggc agtacatcaa tgggcgtgga tagcggtttg actcacgggg	420
	atttccaagt ctccacccca ttgacgtcaa tgggagtttg ttttggcacc aaaatcaacg	480
15	ggactttcca aaatgtcgta acaactccgc cccattgacg caaatgggcg gtaggcgtgt	540
	acggtgggag gtctatataa gcagagctgg tttagtgaac cggcgccagt cctccgattg	600
20	actgagtcgc ccgggtaccc gtgtatccaa taaaccctct tgcagttgca tccgacttgt	660
	ggctctcgctg ttccttggga gggctcctc tgagtgattg actaccgctc agcgggggctc	720
	tttcatttgg gggctcgtcc gggatcgga gaccctgcc cagggaccac cgaccacca	780
25	ccgggaggta agctggccag caacttatct gtgtctgtcc gattgtctag tgtctatgac	840
30		
35		
40		
45		
50		
55		

EP 2 913 402 B1

	tgatTTTTatg	cgcctgcgtc	ggtactagtt	agctaactag	ctctgtatct	ggcggacccg	900
	tggtggaact	gacgagttcg	gaacaccccg	ccgcaaccct	gggagacgtc	ccagggactt	960
5	cgggggcccgt	ttttgtggcc	cgacctgagt	ccaaaaatcc	cgatcgtttt	ggactctttg	1020
	gtgcaccccc	cttagaggag	ggatatgtgg	ttctggtagg	agacgagaac	ctaaaacagt	1080
	tcccgcctcc	gtctgaattt	ttgctttcgg	tttgggaccg	aagccgcgcc	gcgcgctcttg	1140
10	tctgctgcag	catcgttctg	tgttgtctct	gtctgactgt	gtttctgtat	ttgtctgaga	1200
	atatgggcca	gactgttacc	actcccttaa	gtttgacctt	aggtcactgg	aaagatgtcg	1260
	agcggatcgc	tcacaaccag	tcggtagatg	tcaagaagag	acgttgggtt	accttctgct	1320
15	ctgcagaatg	gccaaccttt	aacgtcggat	ggccgcgaga	cggcaccttt	aaccgagacc	1380
	tcatacccca	ggttaagatc	aaggtctttt	cacctggccc	gcatggacac	ccagaccagg	1440
20	tcccctacat	cgtgacctgg	gaagccttgg	cttttgacct	ccctccctgg	gtcaagccct	1500
	ttgtacaccc	taagcctccg	cctcctcttc	ctccatccgc	cccgtctctc	ccccttgaac	1560
	ctcctcgttc	gaccccgctt	cgatcctccc	tttatccagc	cctcactcct	tctctaggcg	1620
25	ccaaacctaa	acctcaagtt	ctttctgaca	gtggggggcc	gctcatcgac	ctacttacag	1680
	aagaccccc	gccttatagg	gacccaagac	caccccttcc	cgacagggac	ggaaatggtg	1740
	gagaagcgac	ccctgcggga	gaggcaccgg	accctcccc	aatggcatct	cgcctacgtg	1800
30	ggagacggga	gccccctgtg	gccgactcca	ctacctcgca	ggcattcccc	ctccgcgcag	1860
	gaggaaacgg	acagcttcaa	tactggccgt	tctcctcttc	tgacctttac	aactggaaaa	1920
35	ataataaccc	ttctttttct	gaagatccag	gtaaaactgac	agctctgatc	gagtctgttc	1980
	tcatacccca	tcagcccacc	tgggacgact	gtcagcagct	gttggggact	ctgctgaccg	2040
	gagaagaaaa	acaacgggtg	ctcttagagg	ctagaaaggc	ggtgcggggc	gatgatgggc	2100
40	gccccactca	actgcccact	gaagtcgatg	ccgcttttcc	cctcgagcgc	ccagactggg	2160
	attacaccac	ccaggcaggt	aggaaccacc	tagtccacta	tcgccagttg	ctcctagcgg	2220
	gtctccaaaa	cgcgggcaga	agccccacca	atttggccaa	ggtaaaagga	ataacacaag	2280
45	ggcccaatga	gtctccctcg	gccttcctag	agagacttaa	ggaagcctat	cgcaggtaca	2340
	ctccttatga	ccctgaggac	ccagggcaag	aaactaatgt	gtctatgtct	ttcatttggc	2400
	agtctgcccc	agacattggg	agaaagttag	agaggttaga	agatttaaaa	aacaagacgc	2460
50	ttggagattt	ggttagagag	gcagaaaaga	tctttaataa	acgagaaacc	ccggaagaaa	2520
	gagaggaacg	tatcaggaga	gaaacagagg	aaaaagaaga	acgccgtagg	acagaggatg	2580
55	agcagaaaga	gaaagaaaga	gatcgtagga	gacatagaga	gatgagcaag	ctatttggcca	2640
	ctgtcgtttag	tggacagaaa	caggatagac	agggaggaga	acgaaggagg	tcccaactcg	2700

EP 2 913 402 B1

	atcgcgacca	gtgtgcctac	tgcaaagaaa	aggggactg	ggctaaagat	tgtcccaaga	2760
	aaccacgagg	acctcgggga	ccaagacccc	agacctccct	cctgacccta	gatgactagg	2820
5	gaggtcaggg	tcaggagccc	ccccctgaac	ccaggataac	cctcaaagtc	ggggggcaac	2880
	ccgtcacctt	cctggtagat	actggggccc	aacactccgt	gctgacccaa	aatcctggac	2940
	ccctaagtga	taagtctgcc	tgggtccaag	gggctactgg	aggaaagcgg	tatcgctgga	3000
10	ccacggatcg	caaagtacat	ctagctaccg	gtaaggtcac	ccactctttc	ctccatgtac	3060
	cagactgtcc	ctatcctctg	ttaggaagag	atttgctgac	taaactaaaa	gcccaaatcc	3120
	actttgaggg	atcaggagcc	caggttatgg	gaccaatggg	gcagcccctg	caagtgttga	3180
15	ccctaaatat	agaagatgag	catcggctac	atgagacctc	aaaagagcca	gatgtttctc	3240
	tagggtccac	atggctgtct	gattttcctc	aggcctgggc	ggaaaccggg	ggcatgggac	3300
	tggcagttcg	ccaagctcct	ctgatcatac	ctctgaaagc	aacctctacc	cccgtgtcca	3360
20	taaaacaata	ccccatgtca	caagaagcca	gactggggat	caagccccac	atacagagac	3420
	tgttggaacca	gggaatactg	gtaccctgcc	agtccccctg	gaacacgccc	ctgctacccg	3480
25	ttaagaaacc	agggactaat	gattataggc	ctgtccagga	tctgagagaa	gtcaacaagc	3540
	gggtggaaga	catccacccc	accgtgccc	acccttacia	cctcttgagc	gggctccac	3600
	cgtcccacca	gtggtacact	gtgcttgatt	taaaggatgc	ctttttctgc	ctgagactcc	3660
30	accccaccag	tcagcctctc	ttcgcccttg	agtggagaga	tccagagatg	ggaatctcag	3720
	gacaattgac	ctggaccaga	ctcccacagg	gtttcaaaaa	cagtcccacc	ctgtttgatg	3780
	aggcactgca	cagagacct	gcagacttcc	ggatccagca	cccagacttg	atcctgctac	3840
35	agtacgtgga	tgacttactg	ctggccgcca	cttctgagct	agactgcca	caaggtactc	3900
	gggccctgtt	acaaacccta	gggaacctcg	ggtatcgggc	ctcggccaag	aaagcccaaa	3960
	tttgccagaa	acaggtcaag	tatctgggg	atcttctaaa	agagggtcag	agatggctga	4020
40	ctgaggccag	aaaagagact	gtgatgggg	agcctactcc	gaagaccctc	cgacaactaa	4080
	gggagttcct	agggacggca	ggcttctgtc	gcctctggat	ccctgggttt	gcagaaatgg	4140
	cagccccctt	gtaccctctc	acaaaaacgg	ggactctgtt	taattggggc	ccagaccaac	4200
45	aaaaggccta	tcaagaaatc	aagcaagctc	ttctaactgc	cccagccctg	gggttgccag	4260
	atttgactaa	gccctttgaa	ctctttgtcg	acgagaagca	gggctacgcc	aaaggtgtcc	4320
50	taacgcaaaa	actgggacct	tggcgctcgg	cggcggccta	cctgtccaaa	aagctagacc	4380
	cagtagcagc	tgggtggccc	ccttgccctac	ggatggtagc	agccattgcc	gtactgacaa	4440
	aggatgcagg	caagctaacc	atgggacagc	cactagtcac	tctggccccc	catgcagtag	4500
55	aggcactagt	caaacaaccc	cccgaccgct	ggctttccaa	cgccccgatg	actcactatc	4560
	aggccttgct	tttgacacg	gaccgggtcc	agttcggacc	ggtggtagcc	ctgaaccggg	4620

EP 2 913 402 B1

	ctacgctgct	cccactgcct	gaggaagggc	tgcaacacaa	ctgccttgat	atcctggccg	4680
	aagcccacgg	aacccgaccc	gacctaacgg	accagccgct	cccagacgcc	gaccacacct	4740
5	ggtacacgga	tggaagcagt	ctcttacaag	agggacagcg	taaggcgga	gctgcggtga	4800
	ccaccgagac	cgaggtaatc	tgggctaaag	ccctgccagc	cgggacatcc	gctcagcggg	4860
	ctgaactgat	agcactcacc	caggccctaa	agatggcaga	aggtaagaag	ctaaatgttt	4920
10	atactgatag	ccgttatgct	tttgctactg	cccatatcca	tggagaaata	tacagaaggc	4980
	gtggggttgct	cacatcagaa	ggcaaagaga	tcaaaaataa	agacgagatc	ttggccctac	5040
	taaaagccct	ctttctgccc	aaaagactta	gcataatcca	ttgtccagga	catcaaaaagg	5100
15	gacacagcgc	cgaggctaga	ggcaaccgga	tggctgacca	agcggcccga	aaggcagcca	5160
	tcacagagac	tccagacacc	tctaccctcc	tcatagaaaa	ttcatcacc	tacacctcag	5220
20	aacattttca	ttacacagtg	actgatataa	aggacctaac	caagttgggg	gccattttatg	5280
	ataaaacaaa	gaagtattgg	gtctaccaag	gaaaacctgt	gatgcctgac	cagtttactt	5340
	ttgaattatt	agactttctt	catcagctga	ctcacctcag	cttctcaaaa	atgaaggctc	5400
25	tcctagagag	aagccacagt	ccctactaca	tgctgaaccg	ggatcgaaca	ctcaaaaata	5460
	tcactgagac	ctgcaaagct	tgtgcacaag	tcaacgccag	caagtctgcc	gttaaacagg	5520
	gaactagggg	ccgcgggcat	cggcccggca	ctcattggga	gatcgatttc	accgagataa	5580
30	agcccggatt	gtatggctat	aaatatcttc	tagtttttat	agataccttt	tctggctgga	5640
	tagaagcctt	cccaaccaag	aaagaaaccg	ccaaggtcgt	aaccaagaag	ctactagagg	5700
35	agatcttccc	caggttcggc	atgcctcagg	tattgggaac	tgacaatggg	cctgccttcg	5760
	tctccaaggt	gagtcagaca	gtggccgcatc	tgttggggat	tgattggaaa	ttacattgtg	5820
	catacagacc	caaagctca	ggccaggtag	aaagaatgaa	tagaaccatc	aaggagactt	5880
40	taactaaatt	aacgcttgca	actggctcta	gagactgggt	gctcctactc	cccttagccc	5940
	tgtaccgagc	ccgcaacacg	ccgggcccc	atggcctcac	cccatatgag	atcttatatg	6000
	gggcaccccc	gccccttgta	aacttccctg	accctgacat	gacaagagtt	actaacagcc	6060
45	cctctctcca	agctcactta	caggctctct	acttagtcca	gcacgaagtc	tggagacctc	6120
	tggcggcagc	ctaccaagaa	caactggacc	gaccggtgg	acctcaccct	taccgagtcg	6180
	gcgacacagt	gtgggtccgc	cgacaccaga	ctaagaacct	agaacctcgc	tggaaaggac	6240
50	cttacacagt	cctgctgacc	acccccaccg	ccctcaaagt	agacggcatc	gcagcttgga	6300
	tacacgccgc	ccacgtgaag	gctgccgacc	ccgggggtgg	accatcctct	agactgacat	6360
55	ggcgcggttca	acgctctcaa	aacccccctca	agataagatt	aacccgtgga	agcccttaat	6420
	agtcatggga	gtcctgttag	gagtagggat	ggcagagagc	ccccatcagg	tctttaatgt	6480

EP 2 913 402 B1

	aacctggaga gtcaccaacc tgatgactgg gcgtaccgcc aatgccacct ccctcctggg	6540
	aactgtacaa gatgccttcc caaaattata ttttgatcta tgtgatctgg tcggagagga	6600
5	gtgggaccct tcagaccagg aaccgtatgt cgggtatggc tgcaagtacc ccgcaggag	6660
	acagcggacc cggacttttg acttttacgt gtgccctggg cataccgtaa agtcgggggtg	6720
	tgggggacca ggagagggct actgtggtaa atgggggtgt gaaaccaccg gacaggctta	6780
10	ctggaagccc acatcatcgt gggacctaata ctcccttaag cgcggtaaca cccctggga	6840
	cacgggatgc tctaaagttg cctgtggccc ctgctacgac ctctccaaag tatccaattc	6900
	cttccaaggg gctactcgag ggggcagatg caaccctcta gtccctagaat tctactgatgc	6960
15	aggaaaaaag gctaactggg acggggcccaa atcgtgggga ctgagactgt accggacagg	7020
	aacagatcct attaccatgt tctccctgac ccggcaggtc cttaatgtgg gaccccgagt	7080
	ccccataggg cccaaccacg tattaccgga ccaaagactc ccttcctcac caatagagat	7140
20	tgtaccggct ccacagccac ctagccccct caataccagt taccctcctt ccactaccag	7200
	tacaccctca acctccccta caagtccaag tgtcccacag ccacccccag gaactggaga	7260
25	tagactacta gctctagtca aaggagccta tcaggcgctt aacctcacca atcccgacaa	7320
	gaccaagaa tgttggtgtg gcttagtgtc gggacctcct tattacgaag gagtagcgg	7380
	cgtgggcact tataccaatc attccaccgc tccggccaac tgtacggcca cttcccaaca	7440
30	taagcttacc ctatctgaag tgacaggaca gggcctatgc atgggggagac tacctaaaac	7500
	tcaccaggcc ttatgtaaca ccacccaaag cgcgggtca ggatcctact accttcgagc	7560
	acccgccgga acaatgtggg cttgcagcac tggattgact ccctgcttgt ccaccacgg	7620
35	gctcaatcta accacagatt attgtgtatt agttgaactc tggcccagag taatttacca	7680
	ctccccgat tatatgtatg gtcagcttga acagcgtaac aaatataaaa gagagccagt	7740
	atcattgacc ctggcccttc tactaggagg attaaccatg ggagggattg cagctggaat	7800
40	agggacgggg accactgcct taattaaaac ccagcagttt gagcagcttc atgccgctat	7860
	ccagacagac ctcaacgaag tcgaaaagtc aattaccaac ctagaaaagt cactgacctc	7920
	gttgtctgaa gtagtcctac agaaccgcag aggcctagat ttgctattcc taaaggagg	7980
45	aggtctctgc gcagccctaa aagaagaatg ttgtttttat gcagaccaca cggggctagt	8040
	gagagacagc atggccaaat taagagaaag gcttaatcag agacaaaaac tatttgagac	8100
	aggccaagga tgggtcgaag ggctgtttta tagatcccc tggtttacca ccttaatctc	8160
50	caccatcatg ggacctctaa tagtactctt actgatctta ctctttggac cttgcattct	8220
	caatcgattg gtccaatttg ttaaagacag gatctcagtg gtccaggctc tggttttgac	8280
55	tcagcaatat caccagctaa aaccataga gtacgagcca tgaacgcgtt actggccgaa	8340
	gccgcttgga ataaggccgg tgtgcgtttg tctatatgtt attttccacc atattgccgt	8400

EP 2 913 402 B1

	cttttgga	caatgtgagggcc	cggaaacctg	gccctgtctt	cttgacgagc	attcctaggg	8460
	gtctttcccc	tctcgccaaa	ggaatgcaag	gtctgttgaa	tgtcgtgaag	gaagcagttc	8520
5	ctctggaagc	ttcttgaaga	caaacaacgt	ctgtagcgac	cctttgcagg	cagcggaacc	8580
	ccccacctgg	cgacaggtgc	ctctgcggcc	aaaagccacg	tgtataagat	acacctgcaa	8640
	aggcggcaca	accccagtgc	cacgttgtga	gttgatagt	tgtggaaaga	gtcaaatggc	8700
10	tctcctcaag	cgtattcaac	aaggggctga	aggatgccca	gaaggtagcc	cattgtatgg	8760
	gatctgatct	ggggcctcgg	tgcacatgct	ttacatgtgt	ttagtcgagg	ttaaaaaac	8820
	gtctaggccc	cccgaaccac	ggggacgtgg	ttttcctttg	aaaaacacga	ttataaatgg	8880
15	tgaccggcgg	catggcctcc	aagtgggatc	aaaagggcat	ggatatcgct	tacgaggagg	8940
	ccctgctggg	ctacaaggag	ggcggcgtgc	ctatcggcgg	ctgtctgatc	aacaacaagg	9000
20	acggcagtgt	gctgggcagg	ggccacaaca	tgaggttcca	gaagggtccc	gccaccctgc	9060
	acggcgagat	ctccaccctg	gagaactgtg	gcaggctgga	gggcaagggtg	tacaaggaca	9120
	ccaccctgta	caccaccctg	tccccttgtg	acatgtgtac	cggcgctatc	atcatgtacg	9180
25	gcatccctag	gtgtgtgatc	ggcgagaacg	tgaacttcaa	gtccaagggc	gagaagtacc	9240
	tgcaaacacg	gggccacgag	gtggtggttg	ttgacgatga	gaggtgtaag	aagctgatga	9300
	agcagttcat	cgacgagagg	cctcaggact	ggttcgagga	tatcggcgag	taagcggccg	9360
30	cagataaaat	aaaagatttt	atttagtctc	cagaaaaagg	ggggaatgaa	agaccccacc	9420
	tgtaggtttg	gcaagctagc	ttaagtaacg	ccattttgca	aggcatggaa	aaatacataa	9480
35	ctgagaatag	agaagttcag	atcaagggtca	ggaacagatg	gaacagctga	atatgggcca	9540
	aacaggatat	ctgtggtaag	cagttcctgc	cccggctcag	ggccaagaac	agatggaaca	9600
	gctgaatatg	ggccaaacag	gatatctgtg	gtaagcagtt	cctgccccgg	ctcagggcca	9660
40	agaacagatg	gtccccagat	gcggtccagc	cctcagcagt	ttctagagaa	ccatcagatg	9720
	tttccagggt	gccccaaagg	cctgaaatga	ccctgtgcct	tatttgaaact	aaccaatcag	9780
	ttcgcttctc	gcttctgttc	gcgcgcttct	gctccccgag	ctcaataaaa	gagcccacaa	9840
45	cccctcactc	ggggcgccag	tcctccgatt	gactgagtcg	cccgggtacc	cgtgtatcca	9900
	ataaaccttc	ttgcagttgc	atccgacttg	tggctctcgct	gttccttggg	agggctctct	9960
	ctgagtgatt	gactaccctg	cagcgggggt	ctttcattac	atgtgagcaa	aaggccagca	10020
50	aaaggccagg	aaccgtaaaa	aggccgcgtt	gctggcgttt	ttccataggc	tccgcccccc	10080
	tgacgagcat	cacaaaaatc	gacgctcaag	tcagagggtg	cgaaacccga	caggactata	10140
55	aagataccag	gcgtttcccc	ctggaagctc	cctcgtgcgc	tctcctgttc	cgaccctgcc	10200
	gcttaccgga	tacctgtccg	cctttctccc	ttcgggaagc	gtggcgcttt	ctcaatgctc	10260

EP 2 913 402 B1

acgctgtagg tatctcagtt cgggtgtaggt cggtcgctcc aagctgggct gtgtgcacga 10320
 accccccggtt cagcccgacc gctgcgcctt atccggtaac tatcgtcttg agtccaaccc 10380
 5 ggtaagacac gacttatcgc cactggcagc agccactggt aacaggatta gcagagcgag 10440
 gtatgtaggc ggtgctacag agttcttgaa gtgggtggcct aactacggct aactagaag 10500
 10 gacagtatctt ggtatctgcg ctctgctgaa gccagttacc ttcggaaaaa gagttggtag 10560
 ctcttgatcc ggcaaaaaa ccaccgctgg tagcgggtgt ttttttgttt gcaagcagca 10620
 gattacgcgc agaaaaaaag gatctcaaga agatcctttg atcttttcta cgggggtctga 10680
 15 cgctcagtgg aacgaaaact cacgttaagg gattttggtc atgagattat caaaaaggat 10740
 cttcacctag atccttttaa attaaaaatg aagtttttaa tcaatctaaa gtatatatga 10800
 gtaaacttgg tctgacagtt accaatgctt aatcagttag gcacctatct cagcgatctg 10860
 20 tctatttcgt tcatccatag ttgcctgact ccccgctcgt tagataacta cgatacggga 10920
 gggcttacca tctggcccca gtgctgcaat gataccgcga gaccacgct caccggctcc 10980
 agatttatca gcaataaacc agccagcccg aagggccgag cgcagaagtg gtcctgcaac 11040
 25 tttatccgcc tccatccagt ctattaattg ttgccgggaa gctagagtaa gtagttcgcc 11100
 agttaatagt ttgcgcaacg ttgttgccat tgctgcaggc atcgtggtgt cacgctcgtc 11160
 gtttggtatg gcttcattca gctccggtt ccaacgatca aggcgagtta catgatcccc 11220
 30 catgttgtgc aaaaaagcgg ttagctcctt cggtcctccg atcgttgtca gaagtaagtt 11280
 ggccgcagtg ttatcactca tggttatggc agcactgcat aattctctta ctgtcatgcc 11340
 35 atccgtaaga tgcttttctg tgactggtga gtactcaacc aagtcattct gagaatagtg 11400
 tatgcggcga ccgagttgct cttgcccggc gtcaacacgg gataataccg cgccacatag 11460
 cagaacttta aaagtgtcga tcattggaac acgttcttcg gggcgaaaac tctcaaggat 11520
 40 cttaccgctg ttgagatcca gttcgatgta acccactcgt gcacccaact gatcttcagc 11580
 atcttttact ttcaccagcg tttctgggtg agcaaaaaca ggaaggcaaa atgccgcaaa 11640
 45 aaagggaata agggcgacac ggaaatgttg aatactcata ctcttccttt ttcaatatta 11700
 ttgaagcatt tatcagggtt attgtctcat gagcggatag atatttgaat gtatttagaa 11760
 aaataaacia ataggggttc cgcgcacatt tccccgaaa gtgccacctg acgtctaaga 11820
 50 aaccattatt atcatgacat taacctataa aaataggcgt atcacgaggc cctttcgtct 11880
 tcaagaattc at 11892

<210> 23

<211> 11892

<212> DNA

<213> Artificial Sequence

EP 2 913 402 B1

<220>

<223> RCR Vector - pAC3-yCD

<400> 23

5

10

15

20

25

30

35

40

45

50

55

EP 2 913 402 B1

	tagttattaa tagtaatcaa ttacgggggtc attagttcat agcccatata tggagttccg	60
	cgttacataa cttacggtaa atggcccgcc tggctgaccg cccaacgacc cccgcccatt	120
5	gacgtcaata atgacgtatg ttcccatagt aacgccaata gggactttcc attgacgtca	180
	atgggtggag tattttacggt aaactgccc cttggcagta catcaagtgt atcatatgcc	240
10	aagtacgccc cctattgacg tcaatgacgg taaatggccc gcctggcatt atgccagta	300
	catgacctta tgggactttc ctacttggca gtacatctac gtattagtca tgcctattac	360
	catggtgatg cggtttttggc agtacatcaa tgggcgtgga tagcggtttg actcacgggg	420
15	atttccaagt ctccacccca ttgacgtcaa tgggagtttg ttttggcacc aaaatcaacg	480
	ggactttcca aaatgtcgta acaactccgc ccattgacg caaatgggcg gtaggcgtgt	540
	acggtgggag gtctatataa gcagagctgg tttagtgaac cggcgccagt cctccgattg	600
20	actgagtcgc ccgggtaccc gtgtatccaa taaaccctct tgcagttgca tccgacttgt	660
	ggtctcgctg ttccttggga gggctcctc tgagtgattg actaccgctc agcgggggtc	720
	tttcatttgg gggctcgtcc gggatcggga gaccctgccc cagggaccac cgaccacca	780
25	ccgggaggta agctggccag caacttatct gtgtctgtcc gattgtctag tgtctatgac	840
	tgattttatg cgcctgcgtc ggtactagtt agctaactag ctctgtatct ggcggacccg	900
30	tgggtggaact gacgagttcg gaacaccccg ccgcaaccct gggagacgtc ccagggactt	960
	cggggggccgt ttttgtggcc cgacctgagt ccaaaaatcc cgatcgtttt ggactctttg	1020
	gtgcaccccc cttagaggag ggatatgtgg ttctggtagg agacgagaac ctaaaacagt	1080
35	tcccgccctcc gtctgaattt ttgctttcgg tttgggaccg aagccgcgcc gcgcgtcttg	1140
	tctgctgcag catcgttctg tgttgtctct gtctgactgt gtttctgtat ttgtctgaga	1200
40	atatgggcca gactgttacc actcccttaa gtttgacctt aggtcactgg aaagatgtcg	1260
	agcggatcgc tcacaaccag tcggtagatg tcaagaagag acgttgggtt accttctgct	1320
	ctgcagaatg gccaaccttt aacgtcggat ggccgcgaga cggcaccttt aaccgagacc	1380
45	tcacacccca ggttaagatc aaggctcttt cacctggccc gcatggacac ccagaccagg	1440
	tcccctacat cgtgacctgg gaagccttgg cttttgaccc ccctccctgg gtcaagccct	1500
	ttgtacaccc taagcctccg cctcctcttc ctccatccgc cccgtctctc ccccttgaac	1560
50	ctcctcgttc gaccccgccct cgatcctccc tttatccagc cctcactcct tctctaggcg	1620
	ccaaacctaa acctcaagtt ctttctgaca gtggggggcc gctcatcgac ctacttacag	1680
	aagaccccc gccttatagg gacccaagac caccctctc cgacagggac ggaaatggtg	1740
55	gagaagcgac ccctgcggga gaggcaccgg acccctcccc aatggcatct cgcctacgtg	1800

EP 2 913 402 B1

	ggagacggga	gccccctgtg	gccgactcca	ctacctcgca	ggcattcccc	ctccgcgcag	1860
	gaggaaacgg	acagcttcaa	tactggccgt	tctcctcttc	tgacctttac	aactggaaaa	1920
5	ataataaccc	ttctttttct	gaagatccag	gtaaactgac	agctctgac	gagtctgttc	1980
	tcatcaccca	tcagcccacc	tgggacgact	gtcagcagct	gttggggact	ctgctgaccg	2040
	gagaagaaaa	acaacgggtg	ctcttagagg	ctagaaaggc	ggtgcggggc	gatgatgggc	2100
10	gccccactca	actgcccaat	gaagtcgatg	ccgcttttcc	cctcgagcgc	ccagactggg	2160
	attacaccac	ccaggcaggt	aggaaccacc	tagtccacta	tcgccagttg	ctcctagcgg	2220
	gtctccaaaa	cgcgggcaga	agccccacca	atttggccaa	ggtaaaagga	ataacacaag	2280
15	ggcccaatga	gtctccctcg	gccttcctag	agagacttaa	ggaagcctat	cgcaggtaca	2340
	ctccttatga	ccctgaggac	ccagggcaag	aaactaatgt	gtctatgtct	ttcatttggc	2400
	agtctgcccc	agacattggg	agaaagttag	agagggttaga	agatttaaaa	aacaagacgc	2460
20	ttggagatth	ggttagagag	gcagaaaaga	tctttaataa	acgagaaacc	ccggaagaaa	2520
	gagaggaacg	tatcaggaga	gaaacagagg	aaaaagaaga	acgccgtagg	acagaggatg	2580
25	agcagaaaga	gaaagaaaga	gatcgtagga	gacatagaga	gatgagcaag	ctattggcca	2640
	ctgtcgttag	tggacagaaa	caggatagac	agggaggaga	acgaaggagg	tcccaactcg	2700
	atcgcgacca	gtgtgcctac	tgcaaagaaa	aggggcactg	ggctaaagat	tgtcccaaga	2760
30	aaccacgagg	acctcgggga	ccaagacccc	agacctccct	cctgacccta	gatgactagg	2820
	gaggtcaggg	tcaggagccc	ccccctgaac	ccaggataac	cctcaaagtc	ggggggcaac	2880
	ccgtcacctt	cctggtagat	actggggccc	aactactccgt	gctgacccaa	aatcctggac	2940
35	ccctaagtga	taagtctgcc	tgggtccaag	gggctactgg	aggaaagcgg	tatcgctgga	3000
	ccacggatcg	caaagtacat	ctagctaccg	gtaaggtcac	ccactctttc	ctccatgtac	3060
	cagactgtcc	ctatcctctg	ttaggaagag	atttgctgac	taaactaaaa	gccc aaatcc	3120
40	actttgaggg	atcaggagcc	caggttatgg	gaccaatggg	gcagcccctg	caagtgttga	3180
	ccctaaatat	agaagatgag	catcggctac	atgagacctc	aaaagagcca	gatgtttctc	3240
	taggggccac	atggctgtct	gattttcctc	aggcctgggc	ggaaaccggg	ggcatgggac	3300
45	tggcagttcg	ccaagctcct	ctgatcatac	ctctgaaagc	aacctctacc	cccgtgtcca	3360
	taaaacaata	ccccatgtca	caagaagcca	gactggggat	caagccccac	atacagagac	3420
50	tgttggaaca	gggaatactg	gtaccctgcc	agtccccctg	gaacacgccc	ctgctacccg	3480
	ttaagaaacc	agggactaat	gattataggc	ctgtccagga	tctgagagaa	gtcaacaagc	3540
	gggtggaaga	catccacccc	accgtgccca	acccttacia	cctcttgagc	gggctccac	3600
55	cgtcccacca	gtggtacact	gtgcttgatt	taaaggatgc	ctttttctgc	ctgagactcc	3660
	accccaccag	tcagcctctc	ttcgcctttg	agtggagaga	tccagagatg	ggaatctcag	3720

EP 2 913 402 B1

	gacaattgac	ctggaccaga	ctcccacagg	gtttcaaaaa	cagtcccacc	ctgtttgatg	3780
	aggcactgca	cagagaccta	gcagacttcc	ggatccagca	cccagacttg	atcctgctac	3840
5	agtacgtgga	tgacttactg	ctggccgcca	cttctgagct	agactgccaa	caaggtactc	3900
	gggccctgtt	acaaacccta	gggaacctcg	ggtatcgggc	ctcggccaag	aaagcccaaa	3960
	tttgccagaa	acaggtcaag	tatctgggggt	atcttctaaa	agagggtcag	agatggctga	4020
10	ctgaggccag	aaaagagact	gtgatggggc	agcctactcc	gaagaccctt	cgacaactaa	4080
	gggagttcct	agggacggca	ggcttctgtc	gcctctggat	ccctgggttt	gcagaaatgg	4140
	cagccccctt	gtaccctctc	acaaaaacgg	ggactctgtt	taattggggc	ccagaccaac	4200
15	aaaaggccta	tcaagaaatc	aagcaagctc	ttctaactgc	cccagccctg	gggttgccag	4260
	atttgactaa	gccctttgaa	ctctttgtcg	acgagaagca	gggctacgcc	aaaggtgtcc	4320
20	taacgcaaaa	actgggacct	tggcgtcggc	cggtggccta	cctgtccaaa	aagctagacc	4380
	cagtagcagc	tgggtggccc	ccttgcctac	ggatggtagc	agccattgcc	gtactgacaa	4440
	aggatgcagg	caagctaacc	atgggacagc	cactagtcat	tctggccccc	catgcagtag	4500
25	aggcactagt	caaacaaccc	cccgaccgct	ggctttccaa	cgcccggatg	actcactatc	4560
	aggccttgct	tttggaacag	gaccgggtcc	agttcggacc	ggtggtagcc	ctgaacccgg	4620
	ctacgctgct	cccactgcct	gaggaagggc	tgcaacacaa	ctgccttgat	atcctggccg	4680
30	aagcccacgg	aacccgaccc	gacctaacgg	accagccgct	cccagacgcc	gaccacacct	4740
	ggtacacgga	tggaagcagt	ctcttacaag	agggacagcg	taaggcggga	gctgcggtga	4800
	ccaccgagac	cgaggtaatc	tgggctaaag	ccctgccagc	cgggacatcc	gctcagcggg	4860
35	ctgaactgat	agcactcacc	caggccctaa	agatggcaga	aggtaagaag	ctaaatgttt	4920
	atactgatag	ccgttatgct	tttgctactg	cccatatcca	tggagaaata	tacagaaggc	4980
40	gtggggttgct	cacatcagaa	ggcaaagaga	tcaaaaataa	agacgagatc	ttggccctac	5040
	taaaagccct	ctttctgccc	aaaagactta	gcataatcca	ttgtccagga	catcaaaaagg	5100
	gacacagcgc	cgaggctaga	ggcaaccgga	tggctgacca	agcggcccga	aaggcagcca	5160
45	tcacagagac	tccagacacc	tctaccctcc	tcatagaaaa	ttcatcacc	tacacctcag	5220
	aacattttca	ttacacagtg	actgatataa	aggacctaac	caagttgggg	gccattttatg	5280
	ataaaacaaa	gaagtattgg	gtctaccaag	gaaaacctgt	gatgcctgac	cagtttactt	5340
50	ttgaattatt	agactttctt	catcagctga	ctcacctcag	cttctcaaaa	atgaaggctc	5400
	tcctagagag	aagccacagt	ccctactaca	tgctgaaccg	ggatcgaaca	ctcaaaaata	5460
55	tcactgagac	ctgcaaagct	tgtgcacaag	tcaacgccag	caagtctgcc	gttaaacagg	5520
	gaactagggt	ccgcgggcat	cggcccggca	ctcattggga	gatcgatttc	accgagataa	5580

EP 2 913 402 B1

	agccccgatt	gtatggctat	aaatatcttc	tagtttttat	agataccttt	tctggctgga	5640
	tagaagcctt	cccaaccaag	aaagaaaccg	ccaaggtcgt	aaccaagaag	ctactagagg	5700
5	agatcttccc	caggttcggc	atgcctcagg	tattgggaac	tgacaatggg	cctgccttcg	5760
	tctccaaggt	gagtcagaca	gtggccgata	tgttggggat	tgattggaaa	ttacattgtg	5820
	catacagacc	ccaaagctca	ggccaggtag	aaagaatgaa	tagaaccatc	aaggagactt	5880
10	taactaaatt	aacgcttgca	actggctcta	gagactgggt	gtcctactc	cccttagccc	5940
	tgtaccgagc	ccgcaacacg	ccgggcccc	atggcctcac	cccatatgag	atcttatatg	6000
	gggcaccccc	gccccttgta	aacttccttg	accctgacat	gacaagagtt	actaacagcc	6060
15	cctctctcca	agctcactta	caggctctct	acttagtcca	gcacgaagtc	tggagacctc	6120
	tggcggcagc	ctaccaagaa	caactggacc	gaccggtggt	acctcaccct	taccgagtcg	6180
	gcgacacagt	gtgggtccgc	cgacaccaga	ctaagaacct	agaacctcgc	tggaaaggac	6240
20	cttacacagt	cctgctgacc	acccccaccg	ccctcaaagt	agacggcatc	gcagcttgga	6300
	tacacgccgc	ccacgtgaag	gctgccgacc	ccgggggtgg	accatcctct	agactgacat	6360
	ggcgcgttca	acgctctcaa	aacccccctca	agataagatt	aaccogtgga	agcccttaat	6420
25	agtcattggga	gtcctgttag	gagtagggat	ggcagagagc	ccccatcagg	tctttaatgt	6480
	aacctggaga	gtcaccaacc	tgatgactgg	gcgtaccgcc	aatgccacct	ccctcctggg	6540
	aactgtacaa	gatgccttcc	caaaattata	ttttgatcta	tgtgatctgg	tcggagagga	6600
30	gtgggaccct	tcagaccagg	aaccgtatgt	cgggtatggc	tgcaagtacc	ccgcagggag	6660
	acagcggacc	cggacttttg	acttttacgt	gtgccctggg	cataccgtaa	agtcgggggtg	6720
35	tgggggacca	ggagagggct	actgtggtaa	atgggggtgt	gaaaccaccg	gacaggctta	6780
	ctggaagccc	acatcatcgt	gggacctaat	ctcccttaag	cgcggttaaca	ccccctggga	6840
	cacgggatgc	tctaaagttg	cctgtggccc	ctgctacgac	ctctccaaag	tatccaattc	6900
40	cttccaaggg	gtactctgag	ggggcagatg	caaccctcta	gtcctagaat	tcactgatgc	6960
	aggaaaaaag	gctaactggg	acgggcccac	atcgtgggga	ctgagactgt	accggacagg	7020
	aacagatcct	attacatgt	tctccctgac	ccggcaggtc	cttaatgtgg	gaccccgagt	7080
45	ccccataggg	cccaaccag	tattaccgca	ccaaagactc	ccttcctcac	caatagagat	7140
	tgtaccggct	ccacagccac	ctagccccct	caataccagt	tacccccctt	ccactaccag	7200
	tacaccctca	acctccccta	caagtccaag	tgtcccacag	ccacccccag	gaactggaga	7260
50	tagactacta	gctctagtca	aaggagccta	tcaggcgctt	aacctcacca	atcccgacaa	7320
	gaccaagaa	tgttggctgt	gcttagtgtc	gggacctcct	tattacgaag	gagtagcggg	7380
	cgtgggcact	tataccaatc	attccaccgc	tccggccaac	tgtacggcca	cttcccaaca	7440
55	taagcttacc	ctatctgaag	tgacaggaca	gggcctatgc	atggggggcag	tacctaaaac	7500

EP 2 913 402 B1

	tcaccaggcc	ttatgtaaca	ccacccaaag	cgccggctca	ggatcctact	accttgacgc	7560
	acccgccgga	acaatgtggg	cttgacgac	tggattgact	ccctgcttgt	ccaccacggt	7620
5	gctcaatcta	accacagatt	attgtgtatt	agttgaactc	tggcccagag	taattttacca	7680
	ctccccgat	tatatgtatg	gtcagcttga	acagcgtacc	aaatataaaa	gagagccagt	7740
	atcattgacc	ctggcccttc	tactaggagg	attaaccatg	ggagggattg	cagctggaat	7800
10	agggacgggg	accactgcct	taattaaaac	ccagcagttt	gagcagcttc	atgccgctat	7860
	ccagacagac	ctcaacgaag	tcgaaaagtc	aattaccaac	ctagaaaagt	cactgacctc	7920
	gttgtctgaa	gtagtcctac	agaaccgcag	aggcctagat	ttgctattcc	taaaggaggg	7980
15	aggtctctgc	gcagccctaa	aagaagaatg	ttgtttttat	gcagaccaca	cggggctagt	8040
	gagagacagc	atggccaaat	taagagaaag	gcttaatcag	agacaaaaac	tattttgagac	8100
20	aggccaagga	tggttcgaag	ggctgtttaa	tagatcccc	tggtttacca	ccttaatctc	8160
	caccatcatg	ggacctctaa	tagtactctt	actgatctta	ctctttggac	cttgcatctc	8220
	caatcgattg	gtccaatttg	ttaaagacag	gatctcagtg	gtccaggctc	tggttttgac	8280
25	tcagcaatat	caccagctaa	aacccataga	gtacgagcca	tgaacgcgtt	actggccgaa	8340
	gccgcttggg	ataaggccgg	tgtgctgttg	tctatatgtt	attttccacc	atattgccgt	8400
	cttttgga	tgtgagggcc	cggaaacctg	gccctgtctt	cttgacgagc	attcctaggg	8460
30	gtctttcccc	tctcgccaaa	ggaatgcaag	gtctgttgaa	tgtcgtgaag	gaagcagttc	8520
	ctctggaagc	ttcttgaaga	caaacaacgt	ctgtagcgac	cctttgcagg	cagcggaacc	8580
35	ccccacctgg	cgacagggtc	ctctgcggcc	aaaagccacg	tgtataagat	acacctgcaa	8640
	aggcggcaca	acccacgtgc	cacgttgtga	gttggtatgt	tgtggaaaga	gtcaaatggc	8700
	tctcctcaag	cgtattcaac	aaggggctga	aggatgccca	gaaggtaccc	cattgtatgg	8760
40	gatctgatct	ggggcctcgg	tgcacatgct	ttacatgtgt	ttagtcgagg	ttaaaaaac	8820
	gtctaggccc	cccgaaccac	ggggacgtgg	ttttcctttg	aaaaacacga	ttataaatgg	8880
	tgacaggggg	aatggcaagc	aagtgggatc	agaagggtat	ggacattgcc	tatgaggagg	8940
45	cggccttagg	ttacaaagag	ggtggtgttc	ctattggcgg	atgtcttatc	aataacaaag	9000
	acggaagtgt	tctcggtcgt	ggtcacaaca	tgagatttca	aaagggatcc	gccacactac	9060
	atggtgagat	ctccactttg	gaaaactgtg	ggagattaga	gggcaaagtg	tacaaagata	9120
50	ccactttgta	tacgacgctg	tctccatgct	acatgtgtac	aggtgccatc	atcatgtatg	9180
	gtattccacg	ctgtgttgct	ggtgagaacg	ttaatttcaa	aagtaagggc	gagaaatatt	9240
55	tacaaactag	aggtcacgag	gttggtgttg	ttgacgatga	gaggtgtaaa	aagatcatga	9300
	aacaatttat	cgatgaaaga	cctcaggatt	ggtttgaaga	tattggtgag	taggcggccg	9360

EP 2 913 402 B1

	cagataaaat	aaaagatttt	atttagtctc	cagaaaaagg	ggggaatgaa	agaccccacc	9420
	tgtaggtttg	gcaagctagc	ttaagtaacg	ccattttgca	aggcatggaa	aaatacataa	9480
5	ctgagaatag	agaagttcag	atcaagggtca	ggaacagatg	gaacagctga	atatgggccca	9540
	aacaggatat	ctgtggttaag	cagttcctgc	cccggctcag	ggccaagaac	agatggaaca	9600
	gctgaatatg	ggccaaacag	gatatctgtg	gtaagcagtt	cctgccccgg	ctcagggcca	9660
10	agaacagatg	gtccccagat	gcggtccagc	cctcagcagt	ttctagagaa	ccatcagatg	9720
	tttccagggg	gccccaaagg	cctgaaatga	ccctgtgcct	tatttgaact	aaccaatcag	9780
	ttcgcttctc	gcttctgttc	gcgcgcttct	gctccccgag	ctcaataaaa	gagcccacaa	9840
15	cccctcactc	ggggcgccag	tcctccgatt	gactgagtcg	cccgggtacc	cgtgtatcca	9900
	ataaaccctc	ttgcagttgc	atccgacttg	tggctctcgt	gttccttggg	aggggtctct	9960
	ctgagtgatt	gactacccgt	cagcgggggt	ctttcattac	atgtgagcaa	aaggccagca	10020
20	aaaggccagg	aaccgtaaaa	aggccgcgtt	gctggcggtt	ttccataggc	tccgcccccc	10080
	tgacgagcat	cacaaaaatc	gacgctcaag	tcagaggtgg	cgaaaccgga	caggactata	10140
	aagataccag	gcgtttcccc	ctggaagctc	cctcgtgcgc	tctcctgttc	cgaccctgcc	10200
25	gcttaccgga	tacctgtccg	cctttctccc	ttcggaagc	gtggcgcttt	ctcaatgctc	10260
	acgctgtagg	tatctcagtt	cggtgtaggt	cgttcgctcc	aagctgggct	gtgtgcacga	10320
	accccccggt	cagcccgacc	gctgcgcctt	atccgtaac	tatcgtcttg	agccaaccc	10380
30	ggtaagacac	gacttatcgc	cactggcagc	agccactggt	aacaggatta	gcagagcgag	10440
	gtatgtaggc	ggtgctacag	agttcttgaa	gtggtggcct	aactacggct	acactagaag	10500
35	gacagtatth	ggtatctgcg	ctctgctgaa	gccagttacc	ttcggaaaaa	gagttggtag	10560
	ctcttgatcc	ggcaaacaaa	ccaccgctgg	tagcgggtgt	ttttttgttt	gcaagcagca	10620
	gattacgcgc	agaaaaaaag	gatctcaaga	agatcctttg	atcttttcta	cggggctctga	10680
40	cgctcagtg	aacgaaaact	cacgttaagg	gattttggtc	atgagattat	caaaaaggat	10740
	cttcacctag	atccttttta	attaaaaatg	aagtttttaa	tcaatctaaa	gtatatatga	10800
	gtaaacttgg	tctgacagtt	accaatgctt	aatcagtgag	gcacctatct	cagcgatctg	10860
45	tctatttcgt	tcatccatag	ttgcctgact	ccccgtcgtg	tagataacta	cgatacggga	10920
	gggcttacca	tctggcccca	gtgctgcaat	gataccgcga	gaccacgct	caccggctcc	10980
	agatttatca	gcaataaacc	agccagccgg	aagggccgag	cgcagaagtg	gtcctgcaac	11040
50	tttatccgcc	tccatccagt	ctattaattg	ttgccgggaa	gctagagtaa	gtagttcgcc	11100
	agttaatagt	ttgcgcaacg	ttgttgccat	tgctgcaggc	atcgtggtgt	cacgctcgtc	11160
	gtttggtatg	gcttcattca	gctccgggtc	ccaacgatca	aggcgagtta	catgatcccc	11220
55	catgttgtgc	aaaaaagcgg	ttagctcctt	cggtcctccg	atcgttgtca	gaagtaagtt	11280

EP 2 913 402 B1

	ggccgcagtg ttatcactca tggttatggc agcactgcat aattctctta ctgtcatgcc	11340
	atccgtaaga tgcttttctg tgactgggtga gtactcaacc aagtcattct gagaatagtg	11400
5	tatgcggcga ccgagttgct cttgcccggc gtcaacacgg gataataccg cgccacatag	11460
	cagaacttta aaagtgctca tcattggaaa acgttcttcg gggcgaaaac tctcaaggat	11520
	cttaccgctg ttgagatcca gttcgatgta acccactcgt gcaccaact gatcttcagc	11580
10	atcttttact ttcaccagcg tttctgggtg agcaaaaaca ggaaggcaaa atgccgcaaa	11640
	aaagggaata agggcgacac ggaaatgttg aatactcata ctcttccttt ttcaatatta	11700
	ttgaagcatt tatcaggggtt attgtctcat gagcggatac atatttgaat gtatttagaa	11760
15	aaataaacia ataggggttc cgcgcacatt tccccgaaaa gtgccacctg acgtctaaga	11820
	aaccattatt atcatgacat taacctataa aaataggcgt atcacgaggc cctttcgtct	11880
20	tcaagaattc at	11892

<210> 24

<211> 12007

<212> DNA

25 <213> Artificial Sequence

<220>

<223> RCR Vector - pACE-CD

30 <400> 24

35

40

45

50

55

EP 2 913 402 B1

	tagttattaa tagtaatcaa ttacggggtc attagttcat agcccatata tggagttccg	60
	cgttacataa cttacggtaa atggcccgcc tggctgaccg cccaacgacc cccgcccatt	120
5	gacgtcaata atgacgtatg ttcccatagt aacgccaata gggactttcc attgacgtca	180
	atgggtggag tatttacggt aaactgccc cttggcagta catcaagtgt atcatatgcc	240
10	aagtacgccc cctattgacg tcaatgacgg taaatggccc gcctggcatt atgccagta	300
	catgacctta tgggactttc ctacttggca gtacatctac gtattagtca tgcctattac	360
	catggtgatg cgggttttggc agtacatcaa tgggcgtgga tagcggtttg actcacgggg	420
15	atttccaagt ctccacccca ttgacgtcaa tgggagtttg ttttggcacc aaaatcaacg	480
	ggactttcca aaatgtcgta acaactccgc cccattgacg caaatgggcg gtaggcgtgt	540
	acggtgggag gtctatataa gcagagctgg tttagtgaac cggcgccagt cctccgattg	600
20	actgagtcgc ccgggtaccc gtgtatccaa taaaccctct tgcagttgca tccgacttgt	660
	ggtctcgctg ttccttggga gggctcctc tgagtgattg actaccgctc agcgggggctc	720
25	tttcatttgg gggctcgtcc gggatcgga gaccctgcc cagggaccac cgaccaccca	780
	ccgggaggta agctggccag caacttatct gtgtctgtcc gattgtctag tgtctatgac	840
30	tgattttatg cgcctgcgtc ggtactagtt agctaactag ctctgtatct ggcggacccg	900
35		
40		
45		
50		
55		

EP 2 913 402 B1

	tggtggaact gacgagttcg gaacaccccg ccgcaaccct gggagacgtc ccagggactt	960
	cgggggccgt ttttgtggcc cgacctgagt ccaaaaatcc cgatcgtttt ggactctttg	1020
5	gtgcaccccc cttagaggag ggatatgtgg ttctggtagg agacgagaac ctaaaacagt	1080
	tcccgctcc gtctgaattt ttgctttcgg tttgggaccg aagccgcgcc gcgcgtcttg	1140
	tctgctgcag catcgttctg tgttgtctct gtctgactgt gtttctgtat ttgtctgaga	1200
10	atatgggcca gactgttacc actcccttaa gtttgacctt aggtcactgg aaagatgtcg	1260
	agcggatcgc tcacaaccag tcggtagatg tcaagaagag acgttggggtt accttctgct	1320
	ctgcagaatg gccaaccttt aacgtcggat ggccgcgaga cggcaccttt aaccgagacc	1380
15	tcatcaccca ggttaagatc aaggtctttt cacctggccc gcatggacac ccagaccagg	1440
	tcccctacat cgtgacctgg gaagccttgg cttttgacct cctccctgg gtcaagccct	1500
	ttgtacaccc taagcctccg cctcctcttc ctccatccgc cccgtctctc ccccttgaac	1560
20	ctcctcgttc gaccccgct cgatcctccc tttatccagc cctcactcct tctctaggcg	1620
	ccaaacctaa acctcaagtt ctttctgaca gtggggggcc gctcatcgac ctacttacag	1680
25	aagaccccc gccttatagg gacccaagac caccctctc cgacagggac ggaaatggtg	1740
	gagaagcgac ccctgcggga gaggcaccgg acccctcccc aatggcatct cgcctacgtg	1800
	ggagacggga gcccctgtg gccgactcca ctacctcgca ggcattcccc ctccgcgcag	1860
30	gaggaaacgg acagcttcaa tactggccgt tctcctcttc tgacctttac aactggaaaa	1920
	ataataaccc ttctttttct gaagatccag gtaaaactgac agctctgac gagtctgttc	1980
	tcatcaccca tcagcccacc tgggacgact gtcagcagct gttggggact ctgctgaccg	2040
35	gagaagaaaa acaacgggtg ctcttagagg ctagaaaggc ggtgcggggc gatgatgggc	2100
	gccccactca actgcccact gaagtcgatg ccgcttttcc cctcgagcgc ccagactggg	2160
	attacaccac ccaggcagg aggaaccacc tagtccacta tcgccagttg ctccctagcg	2220
40	gtctccaaaa cgcgggcaga agccccacca atttggccaa ggtaaaagga ataacacaag	2280
	ggcccaatga gtctccctcg gccttcctag agagacttaa ggaagcctat cgcaggtaca	2340
	ctccttatga ccctgaggac ccagggcaag aaactaatgt gtctatgtct ttcatttggc	2400
45	agtctgcccc agacattggg agaaagttag agaggttaga agatttaaaa aacaagacgc	2460
	ttggagattt ggttagagag gcagaaaaga tctttaataa acgagaaacc ccggaagaaa	2520
50	gagaggaacg tatcaggaga gaaacagagg aaaaagaaga acgccgtagg acagaggatg	2580
	agcagaaaga gaaagaaaga gatcgtagga gacatagaga gatgagcaag ctattggcca	2640
	ctgtcgttag tggacagaaa caggatagac agggaggaga acgaaggagg tcccaactcg	2700
55	atcgcgacca gtgtgcctac tgcaaagaaa aggggcactg ggctaaagat tgtcccaaga	2760
	aaccacgagg acctcgggga ccaagacccc agacctccct cctgacccta gatgactagg	2820

EP 2 913 402 B1

	gaggtcaggg	tcaggagccc	ccccctgaac	ccaggataac	cctcaaagtc	ggggggcaac	2880
	ccgtcacctt	cctggtagat	actggggccc	aacactccgt	gctgacccaa	aatcctggac	2940
5	ccctaagtga	taagtctgcc	tgggtccaag	gggctactgg	aggaaagcgg	tatcgctgga	3000
	ccacggatcg	caaagtacat	ctagctaccg	gtaagggtcac	ccactctttc	ctccatgtac	3060
	cagactgtcc	ctatcctctg	ttaggaagag	atttgctgac	taaactaaaa	gccccaaatcc	3120
10	actttgaggg	atcaggagcc	cagggttatgg	gaccaatggg	gcagcccctg	caagtgttga	3180
	ccctaaatat	agaagatgag	catcggctac	atgagacctc	aaaagagcca	gatgtttctc	3240
	taggggtccac	atggctgtct	gattttcctc	aggcctgggc	ggaaaccggg	ggcatgggac	3300
15	tggcagttcg	ccaagctcct	ctgatcatac	ctctgaaagc	aacctctacc	cccgtgtcca	3360
	taaaacaata	ccccatgtca	caagaagcca	gactggggat	caagccccac	atacagagac	3420
20	tgttggaacca	gggaatactg	gtaccctgcc	agtccccctg	gaacacgccc	ctgctacccg	3480
	ttaagaaacc	agggactaat	gattataggc	ctgtccagga	tctgagagaa	gtcaacaagc	3540
	gggtggaaga	catccacccc	accgtgccc	acccttacia	cctcttgagc	gggctcccac	3600
25	cgtcccacca	gtggtacact	gtgcttgatt	taaaggatgc	ctttttctgc	ctgagactcc	3660
	accaccacag	tcagcctctc	ttcgcccttg	agtggagaga	tccagagatg	ggaatctcag	3720
	gacaattgac	ctggaccaga	ctcccacagg	gtttcaaaaa	cagtcccacc	ctgtttgatg	3780
30	aggcactgca	cagagacct	gcagacttcc	ggatccagca	cccagacttg	atcctgctac	3840
	agtacgtgga	tgacttactg	ctggccgcca	cttctgagct	agactgcca	caagggtactc	3900
	gggccctgtt	acaaacccta	gggaacctcg	ggtatcgggc	ctcggccaag	aaagcccaaa	3960
35	tttgccagaa	acaggtcaag	tatctggggg	atcttctaaa	agaggggtcag	agatggctga	4020
	ctgaggccag	aaaagagact	gtgatggggc	agcctactcc	gaagaccctc	cgacaactaa	4080
40	gggagttcct	agggacggca	ggcttctgtc	gcctctggat	ccctgggttt	gcagaaatgg	4140
	cagccccctt	gtaccctctc	accaaaacgg	ggactctgtt	taattggggc	ccagaccaac	4200
	aaaaggccta	tcaagaaatc	aagcaagctc	ttctaactgc	cccagccctg	gggttgccag	4260
45	atttgactaa	gccctttgaa	ctctttgtcg	acgagaagca	gggctacgcc	aaagggtgtcc	4320
	taacgcaaaa	actgggacct	tggcgtcggc	cggtggccta	cctgtccaaa	aagctagacc	4380
	cagtagcagc	tgggtggccc	ccttgccctac	ggatggtagc	agccattgcc	gtactgacaa	4440
50	aggatgcagg	caagctaacc	atgggacagc	cactagtcac	tctggccccc	catgcagtag	4500
	aggcactagt	caaacaaccc	cccgaccgct	ggctttccaa	cggccggatg	actcactatc	4560
55	aggccttgct	tttggaacag	gaccgggtcc	agttcggacc	ggtggtagcc	ctgaacccgg	4620
	ctacgctgct	cccactgcct	gaggaagggc	tgcaacacaa	ctgccttgat	atcctggccg	4680

EP 2 913 402 B1

	aagcccacgg	aacccgaccc	gacctaacgg	accagccgct	cccagacgcc	gaccacacct	4740
	ggtacacgga	tggaagcagt	ctcttacaag	agggacagcg	taaggcgggga	gctgcggtga	4800
5	ccaccgagac	cgaggtaatc	tgggctaaag	ccctgccagc	cgggacatcc	gctcagcggg	4860
	ctgaactgat	agcactcacc	caggccctaa	agatggcaga	aggtaagaag	ctaaatgttt	4920
	atactgatag	ccgttatgct	tttgctactg	cccatatcca	tggagaaata	tacagaaggc	4980
10	gtgggttgct	cacatcagaa	ggcaaagaga	tcaaaaataa	agacgagatc	ttggccctac	5040
	taaaagccct	ctttctgccc	aaaagactta	gcataatcca	ttgtccagga	catcaaaagg	5100
	gacacagcgc	cgaggctaga	ggcaaccgga	tggctgacca	agcgggccga	aaggcagcca	5160
15	tcacagagac	tccagacacc	tctaccctcc	tcatagaaaa	ttcatcacc	tacacctcag	5220
	aacattttca	ttacacagtg	actgatataa	aggacctaac	caagttgggg	gccatttatg	5280
	ataaaacaaa	gaagtattgg	gtctaccaag	gaaaacctgt	gatgcctgac	cagtttactt	5340
20	ttgaattatt	agactttctt	catcagctga	ctcacctcag	cttctcaaaa	atgaaggctc	5400
	tcctagagag	aagccacagt	ccctactaca	tgctgaaccg	ggatcgaaca	ctcaaaaata	5460
25	tcactgagac	ctgcaaagct	tgtgcacaag	tcaacgccag	caagtctgcc	gttaaacagg	5520
	gaactaggg	ccgcgggcat	cggcccggca	ctcattggga	gatcgatttc	accgagataa	5580
	agcccggtt	gtatggctat	aaatatcttc	tagtttttat	agataccttt	tctggctgga	5640
30	tagaagcctt	cccaaccaag	aaagaaaccg	ccaagtcgt	aaccaagaag	ctactagagg	5700
	agatcttccc	caggttcggc	atgcctcagg	tattgggaac	tgacaatggg	cctgccttcg	5760
	tctccaagg	gagtcagaca	gtggccgac	tggtggggat	tgattggaaa	ttacattgtg	5820
35	catacagacc	ccaaagctca	ggccaggtag	aaagaatgaa	tagaaccatc	aaggagactt	5880
	taactaaatt	aacgcttgca	actggctcta	gagactgggt	gctcctactc	cccttagccc	5940
	tgtaccgagc	ccgcaacacg	ccggggcccc	atggcctcac	cccatatgag	atcttatatg	6000
40	gggcaccccc	gccccttgta	aacttccttg	accctgacat	gacaagagtt	actaacagcc	6060
	cctctctcca	agctcactta	caggctctct	acttagtcca	gcacgaagtc	tggagacctc	6120
	tggcggcagc	ctaccaagaa	caactggacc	gaccggtggt	acctcaccct	taccgagtcg	6180
45	gcgacacagt	gtgggtccgc	cgacaccaga	ctaagaacct	agaacctcgc	tggaaaggac	6240
	cttacacagt	cctgctgacc	acccccaccg	ccctcaaagt	agacggcatc	gcagcttggga	6300
	tacacgccgc	ccacgtgaag	gctgccgacc	ccgggggtgg	accatcctct	agactgacat	6360
50	ggcgcggttca	acgctctcaa	aacccccctca	agataagatt	aaccogtgga	agcccttaat	6420
	agtcatggga	gtcctgttag	gagtagggat	ggcagagagc	ccccatcagg	tctttaatgt	6480
55	aacctggaga	gtcaccaacc	tgatgactgg	gcgtaccgcc	aatgccacct	ccctcctggg	6540
	aactgtacaa	gatgccttcc	caaaattata	ttttgatcta	tgtgatctgg	tcggagagga	6600

EP 2 913 402 B1

	gtgggaccct	tcagaccagg	aaccgtatgt	cgggtatggc	tgcaagtacc	ccgcagggag	6660
	acagcggacc	cggacttttg	acttttacgt	gtgccctggg	cataccgtaa	agtcgggggtg	6720
5	tgggggacca	ggagagggct	actgtggtaa	atgggggtgt	gaaaccaccg	gacaggctta	6780
	ctggaagccc	acatcatcgt	gggacctaat	ctcccttaag	cgcggttaaca	ccccctggga	6840
	cacgggatgc	tctaaagttg	cctgtggccc	ctgctacgac	ctctccaaag	tatccaattc	6900
10	cttccaaggg	gctactcgag	ggggcagatg	caaccctcta	gtcctagaat	tactgatgc	6960
	aggaaaaaag	gctaactggg	acgggcccac	atcgtgggga	ctgagactgt	accggacagg	7020
15	aacagatcct	attaccatgt	tctccctgac	ccggcaggtc	cttaatgtgg	gaccccgagt	7080
	ccccataggg	ccaacccag	tattaccgga	caaagactc	ccttcctcac	caatagagat	7140
	tgtaccggct	ccacagccac	ctagccccct	caataccagt	tacccccctt	ccactaccag	7200
20	tacaccctca	acctccccta	caagtccaag	tgtcccacag	ccacccccag	gaactggaga	7260
	tagactacta	gctctagtca	aaggagccta	tcaggcgctt	aacctcacca	atcccgacaa	7320
	gaccaagaa	tggttggtgt	gcttagtgtc	gggacctcct	tattacgaag	gagtagcggt	7380
25	cgtgggcact	tataccaatc	attccaccgc	tccggccaac	tgtacggcca	cttcccaaca	7440
	taagcttacc	ctatctgaag	tgacaggaca	gggcctatgc	atggggggcag	tacctaaaac	7500
	tcaccaggcc	ttatgtaaca	ccacccaaag	cgccggctca	ggatcctact	accttgacgc	7560
30	acccgccgga	acaatgtggg	cttgacgac	tggattgact	ccctgcttgt	ccaccacggt	7620
	gctcaatcta	accacagatt	attgtgtatt	agttgaactc	tggcccagag	taattttacca	7680
35	ctccccgat	tatatgtatg	gtcagcttga	acagcgtacc	aaatataaaa	gagagccagt	7740
	atcattgacc	ctggcccttc	tactaggagg	attaaccatg	ggagggattg	cagctggaat	7800
	agggacgggg	accactgcct	taattaaaac	ccagcagttt	gagcagcttc	atgccgctat	7860
40	ccagacagac	ctcaacgaag	tcgaaaagtc	aattaccaac	ctagaaaagt	cactgacctc	7920
	gttgtctgaa	gtagtcctac	agaaccgcag	aggcctagat	ttgctattcc	taaaggaggg	7980
	aggtctctgc	gcagccctaa	aagaagaatg	ttgtttttat	gcagaccaca	cggggctagt	8040
45	gagagacagc	atggccaaat	taagagaaag	gcttaatcag	agacaaaaac	tatttgagac	8100
	aggccaagga	tggttcgaag	ggctgtttaa	tagatcccc	tggtttacca	ccttaatctc	8160
	caccatcatg	ggacctctaa	tagtactctt	actgatctta	ctctttggac	cttgcatctt	8220
50	caatcgatta	gtccaatttg	ttaaagacag	gatatcagtg	gtccaggctc	tagttttgac	8280
	tcaacaatat	caccagctga	agcctataga	gtacgagcca	tgacgtacgt	tactggccga	8340
55	agccgcttg	aataaggccg	gtgtgcgttt	gtctatatgt	tattttccac	catattgccg	8400
	tcttttgcca	atgtgagggc	ccggaaacct	ggccctgtct	tcttgacgag	cattcctagg	8460

EP 2 913 402 B1

	ggctctttccc	ctctcgccaa	aggaatgcaa	ggctctgttg	atgtcgtgaa	ggaagcagtt	8520
	cctctggaag	cttcttgaag	acaaacaacg	tctgtagcga	ccctttgcag	gcagcggaac	8580
5	ccccacctg	gcgacaggtg	cctctgcggc	caaaagccac	gtgtataaga	tacacctgca	8640
	aaggcggcac	aaccccagtg	ccacgttgtg	agttggatag	ttgtggaaag	agtcaaattg	8700
	ctctcctcaa	gcgtattcaa	caaggggctg	aaggatgccc	agaaggtacc	ccattgtatg	8760
10	ggatctgatc	tggggcctcg	gtgcacatgc	tttacctgtg	tttagtcgag	gttaaaaaaa	8820
	cgtctaggcc	ccccgaacca	cggggacgtg	gttttccttt	gaaaaacacg	ataataccat	8880
	ggtgacaggg	ggaatggcaa	gcaagtggga	tcagaagggt	atggacattg	cctatgagga	8940
15	ggcggcctta	ggttaciaaag	aggggtgggt	tcctattggc	ggatgtctta	tcaataacaa	9000
	agacggaagt	gttctcggtc	gtggtcacaa	catgagatgt	caaaagggat	ccgccacact	9060
	acatgggtgag	atctccactt	tggaaaactg	tgggagatta	gagggcaaag	tgtacaaaga	9120
20	taccactttg	tatacgacgc	tgtctccatg	cgacatgtgt	acaggtgcc	tcacatcatg	9180
	tggatattcca	cgctgtgttg	tcggtgagaa	cgtaattttc	aaaagtaagg	gcgagaaata	9240
25	tttaciaaact	agaggtcacg	aggttgttgt	tgttgacgat	gagaggtgta	aaaagatcat	9300
	gaaacaattt	atcgatgaaa	gacctcagga	ttggtttgaa	gatattgggt	agtaggcggc	9360
	cgcgccatag	ataaaataaa	agattttatt	tagtctccag	aaaaaggggg	gaatgaaaga	9420
30	ccccacctgt	aggtttggca	agctagctta	agtaacgcc	ttttgcaagg	catggaaaaa	9480
	tacataactg	agaatagaga	agttcagatc	aaggtcagga	acagatggaa	cagctgaata	9540
	tggggccaaac	aggatatctg	tggtaagcag	ttcctgcccc	ggctcagggc	caagaacaga	9600
35	tgggaacagct	gaatatgggc	caaacaggat	atctgtggta	agcagttcct	gccccgggtc	9660
	aggggccaaga	acagatgggt	cccagatgcg	gtccagccct	cagcagtttc	tagagaacca	9720
	tcagatgttt	ccagggtgcc	ccaaggacct	gaaatgacct	tgtgccttgt	ttaaactaac	9780
40	caatcagttc	gcttctcgct	tctgttcgcg	cgcttctgct	ccccgagctc	aataaaaagag	9840
	cccacaaccc	ctcactcggg	gcgccagtc	tccgattgac	tgagtcgccc	gggtacctgt	9900
	gtatccaata	aacctctttg	cagttgcac	cgacttgtgg	tctcgctgtt	ccttgggagg	9960
45	gtctcctctg	agtgattgac	tacctgtcag	cggggtctt	tcatttgggg	gctcgtccgg	10020
	gatcgggaga	cccctgcccc	gggaccaccg	accaccacc	gggaggtaa	ctggctgcct	10080
	cgcgcgtttc	ggtgatgacg	gtgaaaacct	ctgacatgtg	agcaaaaggc	cagcaaaagg	10140
50	ccaggaaccg	taaaaaggcc	gcgttgctgg	cgtttttcca	taggctccgc	ccccctgacg	10200
	agcatcacia	aaatcgacgc	tcaagtcaga	ggtggcgaaa	cccagacagga	ctataaagat	10260
55	accaggcggt	tccccctgga	agctccctcg	tgcgctctcc	tgttccgacc	ctgccgctta	10320
	ccggatacct	gtccgccttt	ctcccttcgg	gaagcgtggc	gctttctcaa	tgctcacgct	10380

EP 2 913 402 B1

	gtaggtatct cagttcgggtg taggtcgttc gctccaagct gggctgtgtg cacgaacccc	10440
	ccgttcagcc cgaccgctgc gccttatccg gtaactatcg tcttgagtcc aacccggtaa	10500
5	gacacgactt atcgccactg gcagcagcca ctggtaacag gattagcaga gcgaggtatg	10560
	taggcgggtgc tacagagttc ttgaagtggg ggcctaacta cggctacact agaaggacag	10620
10	tatttggtat ctgctgctctg ctgaagccag ttaccttcgg aaaaagagtt ggtagctctt	10680
	gatccggcaa acaaaccacc gctggtagcg gtgggtttttt tgtttgcaag cagcagatta	10740
	cgcgagagaaa aaaaggatct caagaagatc ctttgatctt ttctacgggg tctgacgctc	10800
15	agtggaacga aaactcacgt taagggattt tgggtcatgag attatcaaaa aggatcttca	10860
	cctagatcct tttaaattaa aaatgaagtt ttaaataaat ctaaagtata tatgagtaaa	10920
	cttgggtctga cagttaccaaa tgcttaatat gtgaggcacc tatctcagcg atctgtctat	10980
20	ttcgttcatc catagttgcc tgactccccg tcgtgtagat aactacgata cgggagggct	11040
	taccatctgg cccagtgct gcaatgatac cgcgagaccc acgctcaccg gctccagatt	11100
	tatcagcaat aaaccagcca gccggaaggg ccgagcgagc aagtggctct gcaactttat	11160
25	ccgcctccat ccagtctatt aattgttgcc gggaagctag agtaagtagt tcgccagtta	11220
	atagtttgcg caacgttggt gccattgctg caggcatcgt ggtgtcacgc tcgtcgtttg	11280
30	gtatggcttc attcagctcc ggttcccaac gatcaaggcg agttacatga tccccatgt	11340
	tgtgcaaaaa agcgggttagc tccttcggtc ctccgatcgt tgtcagaagt aagttggccg	11400
	cagtgttatc actcatgggt atggcagcac tgcataattc tcttactgtc atgccatccg	11460
35	taagatgctt ttctgtgact ggtgagtact caaccaagtc attctgagaa tagtgtatgc	11520
	ggcgaccgag ttgctcttgc ccggcgtcaa cacgggataa taccgcgcca catagcagaa	11580
	ctttaaaagt gctcatcatt ggaaaacgtt cttcggggcg aaaactctca aggatcttac	11640
40	cgctgttgag atccagttcg atgtaaccca ctctgtcacc caactgatct tcagcatctt	11700
	ttactttcac cagcgtttct ggggtgagcaa aaacaggaag gcaaaatgcc gcaaaaaagg	11760
45	gaataagggc gacacggaaa tgttgaatac tcatactctt cttttttcaa tattattgaa	11820
	gcatttatca gggttattgt ctcatgagcg gatacatatt tgaatgtatt tagaaaaata	11880
	aacaaatagg ggttccgcgc acatttcccc gaaaagtgcc acctgacgtc taagaaacca	11940
50	ttattatcat gacattaacc tataaaaaata ggcgtatcac gaggcccttt cgtcttcaag	12000
	aattcat	12007
55	<210> 25	
	<211> 11893	
	<212> DNA	
	<213> Artificial Sequence	

EP 2 913 402 B1

<220>

<223> RCR Vector - pAC3-yCD2

<400> 25

5

10

15

20

25

30

35

40

45

50

55

EP 2 913 402 B1

	tagttattaa tagtaatcaa ttacgggggtc attagttcat agcccatata tggagttccg	60
	cgttacataa cttacggtaa atggcccgcc tggctgaccg cccaacgacc cccgcccatt	120
5	gacgtcaata atgacgtatg ttcccatagt aacgccaata gggactttcc attgacgtca	180
	atgggtggag tatttacggt aaactgccc cttggcagta catcaagtgt atcatatgcc	240
	aagtacgccc cctattgacg tcaatgacgg taaatggccc gcctggcatt atgccagta	300
10	catgacctta tgggactttc ctacttggca gtacatctac gtattagtca tgcctattac	360
	catggtgatg cgggttttggc agtacatcaa tgggcgtgga tagcggtttg actcacgggg	420
15	atttccaagt ctccacccca ttgacgtcaa tgggagtttg ttttggcacc aaaatcaacg	480
	ggactttcca aaatgtcgt acaactccgc ccattgacg caaatgggcg gtaggcgtgt	540
	acggtgggag gtctatataa gcagagctgg tttagtgaac cggcgccagt cctccgattg	600
20	actgagtcgc ccgggtaccc gtgtatccaa taaaccctct tgcagttgca tccgacttgt	660
	ggtctcgctg ttccttggga gggctcctc tgagtgattg actaccgctc agcgggggtc	720
	tttcatttgg gggctcgtcc gggatcggga gaccctgcc cagggaccac cgaccacca	780
25	ccgggaggta agctggccag caacttatct gtgtctgtcc gattgtctag tgtctatgac	840
	tgattttatg cgcctgcgtc ggtactagtt agctaactag ctctgtatct ggcggacccg	900
30	tggtggaact gacgagttcg gaacaccccg ccgcaaccct gggagacgtc ccagggactt	960
	cggggggcgt ttttgtggcc cgacctgagt ccaaaaatcc cgatcgtttt ggactctttg	1020
	gtgcaccccc cttagaggag ggatatgtgg ttctggtagg agacgagaac ctaaaacagt	1080
35	tcccgccctc gtctgaattt ttgctttcgg tttgggaccg aagccgcgcc gcgcgtcttg	1140
	tctgctgcag catcgttctg tgttgtctct gtctgactgt gtttctgtat ttgtctgaaa	1200
	atatgggcca gactgttacc actcccttaa gtttgacctt aggtcactgg aaagatgtcg	1260
40	agcggatcgc tcacaaccag tcggtagatg tcaagaagag acgttgggtt accttctgct	1320
	ctgcagaatg gccaaccttt aacgtcggat ggccgcgaga cggcaccttt aaccgagacc	1380
45	tcacaccca ggttaagatc aaggtctttt cacctggccc gcatggacac ccagaccagg	1440
	tcccctacat cgtgacctgg gaagccttgg cttttgacct ccctccctgg gtcaagccct	1500
	ttgtacaccc taagcctccg cctcctcttc ctccatccgc cccgtctctc ccccttgaac	1560
50	ctcctcgttc gaccccgcc cgcctcctcc tttatccagc cctcactcct tctctaggcg	1620
	ccaaacctaa acctcaagtt ctttctgaca gtggggggcc gctcatcgac ctacttacag	1680
	aagaccccc gccttatagg gacccaagac caccctctc cgacagggac ggaaatggtg	1740
55	gagaagcgac ccctgcggga gaggcaccgg acccctcccc aatggcatct cgcctacgtg	1800

EP 2 913 402 B1

	ggagacggga	gccccctgtg	gccgactcca	ctacctcgca	ggcattcccc	ctccgcgcag	1860
	gaggaaacgg	acagcttcaa	tactggccgt	tctcctcttc	tgacctttac	aactggaaaa	1920
5	ataataaccc	ttctttttct	gaagatccag	gtaaaactgac	agctctgatc	gagtctgtcc	1980
	tcatcaccca	tcagcccacc	tgggacgact	gtcagcagct	gttggggact	ctgctgaccg	2040
	gagaagaaaa	acaacgggtg	ctcttagagg	ctagaaaggc	ggtgcggggc	gatgatgggc	2100
10	gccccactca	actgcccact	gaagtcgatg	ccgcttttcc	cctcgagcgc	ccagactggg	2160
	attacaccac	ccaggcaggt	aggaaccacc	tagtccacta	tcgccagttg	ctcctagcgg	2220
	gtctccaaaa	cgcgggcaga	agccccacca	atttggccaa	ggtaaaagga	ataacacaag	2280
15	ggcccaatga	gtctccctcg	gccttcctag	agagacttaa	ggaagcctat	cgcaggtaca	2340
	ctccttatga	ccctgaggac	ccagggcaag	aaactaatgt	gtctatgtct	ttcatttggc	2400
20	agtctgcccc	agacattggg	agaaagttag	agaggtaga	agatttaaaa	aacaagacgc	2460
	ttggagattt	ggttagagag	gcagaaaaga	tctttaataa	acgagaaacc	ccggaagaaa	2520
	gagaggaacg	tatcaggaga	gaaacagagg	aaaaagaaga	acgccgtagg	acagaggatg	2580
25	agcagaaaga	gaaagaaaga	gatcgtagga	gacatagaga	gatgagcaag	ctattggcca	2640
	ctgtcgtagg	tggacagaaa	caggatagac	agggaggaga	acgaaggagg	tcccaactcg	2700
	atcgcgacca	gtgtgcctac	tgcaaagaaa	aggggcactg	ggctaaagat	tgtcccaaga	2760
30	aaccacgagg	acctcgggga	ccaagacccc	agacctccct	cctgacccta	gatgactagg	2820
	gaggtcaggg	tcaggagccc	ccccctgaac	ccaggataac	cctcaaagtc	ggggggcaac	2880
35	ccgtcacctt	cctggtagat	actggggccc	aacactccgt	gctgacccaa	aatcctggac	2940
	ccctaagtga	taagtctgcc	tgggtccaag	gggctactgg	aggaaagcgg	tatcgctgga	3000
	ccacggatcg	caaagtacat	ctagctaccg	gtaaggtcac	ccactctttc	ctccatgtac	3060
40	cagactgtcc	ctatcctctg	ttaggaagag	atttgctgac	taaactaaaa	gccccaatcc	3120
	actttgaggg	atcaggagcc	caggttatgg	gaccaatggg	gcagcccctg	caagtgttga	3180
	ccctaaatat	agaagatgag	tatcggctac	atgagacctc	aaaagagcca	gatgtttctc	3240
45	taggggtccac	atggctgtct	gattttcctc	aggcctgggc	ggaaaccggg	ggcatgggac	3300
	tggcagttcg	ccaagctcct	ctgatcatac	ctctgaaagc	aacctctacc	cccgtgtcca	3360
	taaaacaata	ccccatgtca	caagaagcca	gactggggat	caagccccac	atacagagac	3420
50	tgttggaacca	gggaatactg	gtaccctgcc	agtccccctg	gaacacgccc	ctgctacccg	3480
	ttaagaaacc	agggactaat	gattataggc	ctgtccagga	tctgagagaa	gtcaacaagc	3540
55	gggtggaaga	catccacccc	accgtgccca	acccttacia	cctcttgagc	gggtccccac	3600
	cgtcccacca	gtggtacact	gtgcttgatt	taaaggatgc	ctttttctgc	ctgagactcc	3660

EP 2 913 402 B1

	accccaccag	tcagcctctc	ttcgccctttg	agtggagaga	tccagagatg	ggaatctcag	3720
	gacaattgac	ctggaccaga	ctcccacagg	gtttcaaaaa	cagtcccacc	ctgtttgatg	3780
5	aggcactgca	cagagaccta	gcagacttcc	ggatccagca	cccagacttg	atcctgctac	3840
	agtacgtgga	tgacttactg	ctggccgcca	cttctgagct	agactgccaa	caaggtactc	3900
	gggccctggt	acaaacccta	gggaacctcg	ggtatcgggc	ctcggccaag	aaagcccaaa	3960
10	tttgccagaa	acaggtcaag	tatctggggg	atcttctaaa	agaggggtcag	agatggctga	4020
	ctgaggccag	aaaagagact	gtgatggggc	agcctactcc	gaagaccctc	cgacaactaa	4080
	gggagttcct	agggacggca	ggcttctgtc	gcctctggat	ccctgggttt	gcagaaatgg	4140
15	cagccccctt	gtaccctctc	acaaaaacgg	ggactctgtt	taattggggc	ccagaccaac	4200
	aaaaggccta	tcaagaaatc	aagcaagctc	ttctaactgc	cccagccctg	gggttgccag	4260
	atttgactaa	gccctttgaa	ctctttgtcg	acgagaagca	gggctacgcc	aaaggtgtcc	4320
20	taacgcaaaa	actgggacct	tggcgctcgg	cgggtggccta	cctgtccaaa	aagctagacc	4380
	cagtagcagc	tgggtggccc	ccttgcctac	ggatggtagc	agccattgcc	gtactgacaa	4440
25	aggatgcagg	caagctaacc	atgggacagc	cactagtcac	tctggccccc	catgcagtag	4500
	aggcactagt	caaacaaccc	cccgaccgct	ggctttccaa	cgcccggatg	actcactatc	4560
	aggccttgct	tttgggacacg	gaccgggtcc	agttcggacc	ggtggtagcc	ctgaaccggg	4620
30	ctacgctgct	cccactgcct	gaggaagggc	tgcaacacaa	ctgccttgat	atcctggccg	4680
	aagcccacgg	aacccgaccc	gacctaacgg	accagccgct	cccagacgcc	gaccacacct	4740
	ggtacacgga	tggaaagcagt	ctcttacaag	agggacagcg	taaggcggga	gctgcggtga	4800
35	ccaccgagac	cgaggtaatc	tgggctaaag	ccctgccagc	cgggacatcc	gctcagcggg	4860
	ctgaactgat	agcactcacc	caggccctaa	agatggcaga	aggtaagaag	ctaaatgttt	4920
	atactgatag	ccgttatgct	tttgctactg	cccatatcca	tggagaaata	tacagaaggc	4980
40	gtgggttgct	cacatcagaa	ggcaaagaga	tcaaaaataa	agacgagatc	ttggccctac	5040
	taaaagccct	ctttctgccc	aaaagactta	gcataatcca	ttgtccagga	catcaaaagg	5100
	gacacagcgc	cgaggctaga	ggcaaccgga	tggctgacca	agcggcccca	aaggcagcca	5160
45	tcacagagac	tccagacacc	tctaccctcc	tcatagaaaa	ttcatcacc	tacacctcag	5220
	aacattttca	ttacacagtg	actgatataa	aggacctaac	caagttgggg	gccatttatg	5280
	ataaaacaaa	gaagtattgg	gtctaccaag	gaaaacctgt	gatgcctgac	cagtttactt	5340
50	ttgaattatt	agactttctt	catcagctga	ctcacctcag	cttctcaaaa	atgaaggctc	5400
	tcctagagag	aagccacagt	ccctactaca	tgctgaaccg	ggatcgaaca	ctcaaaaata	5460
55	tcactgagac	ctgcaaagct	tgtgcacaag	tcaacgccag	caagtctgcc	gttaaacagg	5520
	gaactagggt	ccgcgggcat	cggcccgcca	ctcattggga	gatcgatttc	accgagataa	5580

EP 2 913 402 B1

	agcccgatt	gtatggctat	aaatatcttc	tagtttttat	agataccttt	tctggctgga	5640
	tagaagcctt	cccaaccaag	aaagaaaccg	ccaaggtcgt	aaccaagaag	ctactagagg	5700
5	agatcttccc	caggttcggc	atgcctcagg	tattgggaac	tgacaatggg	cctgccttcg	5760
	tctccaaggt	gagtcagaca	gtggccgatc	tgttggggat	tgattggaaa	ttacattgtg	5820
	catacagacc	caaagctca	ggccaggtag	aaagaatgaa	tagaaccatc	aaggagactt	5880
10	taactaaatt	aacgcttgca	actggctcta	gagactgggt	gctcctactc	cccttagccc	5940
	tgtaccgagc	ccgcaacacg	ccgggcccc	atggcctcac	cccatatgag	atcttatatg	6000
	gggcaccccc	gcccttgta	aacttccttg	accctgacat	gacaagagtt	actaacagcc	6060
15	cctctctcca	agctcactta	caggctctct	acttagtcca	gcacgaagtc	tggagacctc	6120
	tggcggcagc	ctaccaagaa	caactggacc	gaccggtggt	acctcaccct	taccgagtcg	6180
20	gcgacacagt	gtgggtccgc	cgacaccaga	ctaagaacct	agaacctcgc	tggaaaggac	6240
	cttacacagt	cctgctgacc	acccccaccg	ccctcaaagt	agacggcatc	gcagcttggg	6300
	tacacgccgc	ccacgtgaag	gctgccgacc	ccgggggtgg	accatcctct	agactgacat	6360
25	ggcgcgttca	acgctctcaa	aacccccctca	agataagatt	aacccgtgga	agcccttaat	6420
	agtcattgga	gtcctgttag	gagtagggat	ggcagagagc	ccccatcagg	tctttaatgt	6480
	aacctggaga	gtcaccaacc	tgatgactgg	gcgtaccgcc	aatgccacct	ccctcctggg	6540
30	aactgtacaa	gatgccttcc	caaaattata	ttttgatcta	tgtgatctgg	tcggagagga	6600
	gtgggaccct	tcagaccagg	aaccgtatgt	cgggtatggc	tgcaagtacc	ccgcagggag	6660
35	acagcggacc	cggacttttg	acttttacgt	gtgccctggg	cataccgtaa	agtcgggggtg	6720
	tgggggacca	ggagagggct	actgtggtaa	atgggggtgt	gaaaccaccg	gacagggctta	6780
	ctggaagccc	acatcatcgt	gggacctaat	ctcccttaag	cgcggttaaca	ccccctggga	6840
40	cacgggatgc	tctaaagttg	cctgtggccc	ctgctacgac	ctctccaaag	tatccaattc	6900
	cttccaaggg	gctactcgag	ggggcagatg	caaccctcta	gtcctagaat	tactgatgc	6960
	aggaaaaaag	gctaactggg	acgggcccaa	atcgtgggga	ctgagactgt	accggacagg	7020
45	aacagatcct	attaccatgt	tctccctgac	ccggcaggtc	cttaatgtgg	gaccccgagt	7080
	ccccataggg	cccaaccag	tattaccgga	ccaaagactc	ccttcctcac	caatagagat	7140
50	tgtaccggct	ccacagccac	ctagccccct	caataccagt	tacccccctt	ccactaccag	7200
	tacaccctca	acctccccta	caagtccaag	tgtcccacag	ccacccccag	gaactggaga	7260
	tagactacta	gctctagtca	aaggagccta	tcaggcgctt	aacctacca	atcccgacaa	7320
55	gaccaagaa	tgttggctgt	gcttagtgct	gggacctcct	tattacgaag	gagtagcggg	7380
	cgtgggcact	tataccaatc	attccaccgc	tccggccaac	tgtacggcca	cttccaaca	7440

EP 2 913 402 B1

	taagcttacc ctatctgaag tgacaggaca gggcctatgc atgggggagac tacctaaaac	7500
	tcaccaggcc ttatgtaaca ccacccaaag cgccggctca ggatcctact accttgacgc	7560
5	acccgccgga acaatgtggg cttgcagcac tggattgact ccctgcttgt ccaccacggt	7620
	gctcaatcta accacagatt attgtgtatt agttgaactc tggcccagag taattttacca	7680
	ctccccgat tatatgtatg gtcagcttga acagcgtaac aaatataaaa gagagccagt	7740
10	atcattgacc ctggcccttc tactaggagg attaaccatg ggagggattg cagctggaat	7800
	agggacgggg accactgcct taattaaaac ccagcagttt gagcagcttc atgccgctat	7860
	ccagacagac ctcaacgaag tcgaaaagtc aattaccaac ctagaaaagt cactgacctc	7920
15	gttgtctgaa gtagtcctac agaaccgcag aggcctagat ttgctattcc taaaggaggg	7980
	aggtctctgc gcagccctaa aagaagaatg ttgtttttat gcagaccaca cggggctagt	8040
	gagagacagc atggccaaat taagagaaag gcttaatcag agacaaaaac tatttgagac	8100
20	aggccaagga tggttcgaag ggctgtttta tagatcccc tggtttacca ccttaatctc	8160
	caccatcatg ggacctctaa tagtactctt actgatctta ctctttggac cttgcattct	8220
25	caatcgattg gtccaatttg ttaaagacag gatctcagtg gtccaggctc tggttttgac	8280
	tcagcaatat caccagctaa aaccataga gtacgagcca tgaacgcgtt actggccgaa	8340
	gccgcttgga ataaggccgg tgtgcgtttg tctatatgtt attttccacc atattgccgt	8400
30	cttttggaac tgtgagggcc cggaacctg gccctgtctt cttgacgagc attcctaggg	8460
	gtctttcccc tctcgccaaa ggaatgcaag gtctgttgaa tgtcgtgaag gaagcagttc	8520
	ctctggaagc ttcttgaaaga caaacaacgt ctgtagcgac cctttgcagg cagcggaacc	8580
35	ccccacctgg cgacaggtgc ctctgcggcc aaaagccacg tgtataagat acacctgcaa	8640
	aggcggcaca accccagtgc cacgttgtga gttggatagt tgtggaaaga gtcaaatggc	8700
	tctcctcaag cgtattcaac aaggggctga aggatgccca gaaggtaacc cattgtatgg	8760
40	gatctgatct ggggcctcgg tgcacatgct ttacatgtgt ttagtcgagg ttaaaaaaac	8820
	gtctaggccc cccgaaccac ggggacgtgg ttttcctttg aaaaacacga ttataaatgg	8880
	tgaccggcgg catggcctcc aagtgggatc aaaagggcat ggatatcgct tacgaggagg	8940
45	ccctgctggg ctacaaggag ggcggcgtgc ctatcggcgg ctgtctgatc aacaacaagg	9000
	acggcagtgt gctgggcagg ggccacaaca tgagggtcca gaagggtcc gccaccctgc	9060
	acggcgagat ctccaccctg gagaactgtg gcaggctgga gggcaagggtg tacaaggaca	9120
50	ccaccctgta caccaccctg tccccttggtg acatgtgtac cggcgctatc atcatgtacg	9180
	gcatccctag gtgtgtgatc ggcgagaacg tgaacttcaa gtccaagggc gagaagtacc	9240
55	tgcaaaccag gggccacgag gtggtgggtt ttgacgatga gaggtgtaag aagctgatga	9300
	agcagttcat cgacgagagg cctcaggact ggttcgagga tatcggcgag taagcggccg	9360

EP 2 913 402 B1

	cagataaaat	aaaagatttt	atttagtctc	cagaaaaagg	ggggaatgaa	agaccccacc	9420
	tgtaggtttg	gcaagctagc	ttaagtaacg	ccattttgca	aggcatggaa	aaatacataa	9480
5	ctgagaatag	agaagttcag	atcaaggtca	ggaacagatg	gaacagctga	atatgggcca	9540
	aacaggatat	ctgtggtaag	cagttcctgc	cccggctcag	ggccaagaac	agatggaaca	9600
	gctgaatatg	ggccaaacag	gatatctgtg	gtaagcagtt	cctgccccgg	ctcagggcca	9660
10	agaacagatg	gtccccagat	gcggtccagc	cctcagcagt	ttctagagaa	ccatcagatg	9720
	tttccagggt	gccccaaagg	cctgaaatga	ccctgtgcct	tatttgaact	aaccaatcag	9780
15	ttcgcttctc	gcttctgttc	gcgcgcttct	gctccccgag	ctcaataaaa	gagcccacaa	9840
	cccctcactc	ggggcgccag	tcctccgatt	gactgagtcg	cccgggtacc	cgtgtatcca	9900
	ataaaccctc	ttgcagttgc	atccgacttg	tggctctcgt	gttccttggg	agggctctct	9960
20	ctgagtgatt	gactaccctg	cagcgggggt	ctttcattac	atgtgagcaa	aaggccagca	10020
	aaaggccagg	aaccgtaaaa	aggccgcgtt	gctggcgttt	ttccataggc	tccgcccccc	10080
	tgacgagcat	cacaaaaatc	gacgctcaag	tcagaggtgg	cgaaacccga	caggactata	10140
25	aagataccag	gcgtttcccc	ctggaagctc	cctcgtgcgc	tctcctgttc	cgaccctgcc	10200
	gcttaccgga	tacctgtccg	cctttctccc	ttcgggaagc	gtggcgcttt	ctcatagctc	10260
	acgctgtagg	tatctcagtt	cggtgtaggt	cgttcgctcc	aagctgggct	gtgtgcacga	10320
30	acccccggt	cagcccgacc	gctgcgcctt	atccggtaac	tatcgtcttg	agtccaaccc	10380
	ggtaagacac	gacttatcgc	cactggcagc	agccactggt	aacaggatta	gcagagcgag	10440
35	gtatgtaggc	ggtgctacag	agttcttgaa	gtggtggcct	aactacggct	acactagaag	10500
	gacagtatth	ggtatctgcg	ctctgctgaa	gccagttacc	ttcggaaaaa	gagttggtag	10560
	ctcttgatcc	ggcaaacaaa	ccaccgctgg	tagcgggtgt	ttttttgttt	gcaagcagca	10620
40	gattacgcgc	agaaaaaaag	gatctcaaga	agatcctttg	atcttttcta	cggggtctga	10680
	cgctcagtg	aacgaaaact	cacgttaagg	gattttggtc	atgagattat	caaaaaggat	10740
	cttcacctag	atccttttaa	attaaaaatg	aagttttaaa	tcaatctaaa	gtatatatga	10800
45	gtaaacttgg	tctgacagtt	accaatgctt	aatcagtgag	gcacctatct	cagcgatctg	10860
	tctatttcgt	tcattccatag	ttgcctgact	ccccgtcgtg	tagataacta	cgatacggga	10920
50	gggcttacca	tctggcccca	gtgctgcaat	gataccgcga	gaccacgct	caccggctcc	10980
	agatttatca	gcaataaacc	agccagccgg	aagggccgag	cgcagaagtg	gtcctgcaac	11040
	tttatccgcc	tccatccagt	ctattaattg	ttgccgggaa	gctagagtaa	gtagttcgcc	11100
55	agttaatatg	ttgcgcaacg	ttgttgccat	tgctgcaggc	atcgtggtgt	cacgctcgtc	11160
	gtttggtatg	gcttcattca	gctccggttc	ccaacgatca	aggcgagtta	catgatcccc	11220

	catgttgtgc aaaaaagcgg ttagctcctt cggtcctccg atcgttgtca gaagtaagtt	11280
	ggccgcagtg ttatcactca tggttatggc agcactgcat aattctctta ctgtcatgcc	11340
5	atccgtaaga tgcttttctg tgactggtga gtactcaacc aagtcattct gagaatagtg	11400
	tatgcggcga ccgagttgct cttgcccggc gtcaacacgg gataataccg cgccacatag	11460
	cagaacttta aaagtgctca tcattggaaa acgttcttcg gggcgaaaac tctcaaggat	11520
10	cttaccgctg ttgagatcca gttcgatgta acccactcgt gcacccaact gatcttcagc	11580
	atcttttact ttcaccagcg tttctgggtg agcaaaaaca ggaaggcaaa atgccgcaaa	11640
	aaaggggaata agggcgacac ggaaatgttg aatactcata ctcttccttt ttcaatatta	11700
15	ttgaagcatt tatcaggggtt attgtctcat gagcggatac atatttgaat gtatttagaa	11760
	aaataaacia ataggggttc cgcgcacatt tccccgaaaa gtgccacctg acgtctaaga	11820
20	aaccattatt atcatgacat taacctataa aaataggcgt atcacgaggc cctttcgtct	11880
	tcaagaattc cat	11893

25 Claims

1. A composition comprising a primer pair and probe selected from the group consisting of:

- (a) SEQ ID NO: 7, 8 and 9; and
 (b) SEQ ID NO: 10, 11 and 12.

2. A method of determining viral content in an *in vitro* sample from a subject about to undergo or undergoing a retroviral gene delivery therapy using an MLV-related virus, comprising:

- contacting the sample with one or more primer pairs as set forth in claim 1 under conditions suitable for nucleic acid amplification to obtain amplified products;
 contacting the sample with a one or more probes that hybridizes to the amplified product;
 detecting a hybridized product;
 indicating that the subject has viral content comprising an MLV-related virus.

3. The method of claim 2, wherein the MLV-related virus is a recombinant retroviral vector used in gene therapy, and/or wherein the method is performed prior to a therapeutic regimen comprising a gene therapy vector, and/or wherein the method is carried out prior to delivery of a MLV-related retroviral vector for gene delivery, and/or wherein the method is carried out following delivery of a MLV-related retroviral vector for gene delivery, and/or wherein the MLV-related virus comprises a 5' LTR, gag, pol, env genes, a regulatory domain 3' of the env gene linked to a heterologous polynucleotide to be delivered and a 3' LTR and a promoter for expression in mammalian cells in the 5'-LTR, preferably wherein the regulatory domain is an internal ribosome entry site (IRES), and/or preferably wherein the heterologous polynucleotide encodes a polypeptide having cytosine deaminase activity,
 and/or wherein the method monitors the spread of the MLV-related retroviral vector, preferably wherein the method is carried out routinely over the course years.

4. A method for detecting the presence of a viral agent in a sample comprising:

- measuring the amount of a polynucleotide in a sample using a quantitative polymerase chain reaction or other amplification process comprising a oligonucleotide primer pair and probe selected from:

- (i) SEQ ID NO: 7, 8 and 9; and

(ii) SEQ ID NO: 10, 11 and 12.

5. The method of claim 4, wherein the polynucleotide is DNA,
and/or wherein the polynucleotide is RNA,
and/or wherein the quantitative polymerase chain reaction is RT-qPCR,
and/or wherein measuring detects a single copy of a viral agent related nucleic acid,
and/or wherein the viral agent comprises a MLV-related virus,
and/or wherein the sample is mammalian tissue,
and/or wherein the sample is mammalian blood,
and/or wherein the viral agent is a gene therapy vector, preferably wherein the gene therapy vector is a replication-competent vector,
and/or wherein the method is performed prior to a therapeutic regimen comprising a gene therapy vector treatment,
and/or wherein the method is performed to monitor the dosage of a therapeutic regimen comprising a gene therapy vector in a subject,
and/or wherein the gene therapy vector comprises a replication competent MLV vector,
and/or wherein the method is performed subsequent to a therapeutic regimen comprising a gene therapy vector,
and/or wherein the method is performed to monitor the dosage of a therapeutic regimen comprising a gene therapy vector,
and/or wherein the sample is from a human subject.
6. A kit for performing the method of any one of claims 2 to 5, the kit comprising a primer pair as defined in claim 1.
7. A method of any one of claims 2 to 5 for detecting <100 copies of MLV-related DNA in a sample extracted from fixed histopathological sections,
and/or for detecting <100 copies of MLV-related RNA in a sample extracted from fixed histopathological sections.
8. A method of any of claims 2 to 5, wherein the method is a method for detecting MLV-related viruses in plasma or serum from a mammalian host.
9. A method of screening a blood supply or tissue bank for infection by a MLV-related virus comprising an amplification reaction on the blood supply or tissue bank utilizing a primer pair as set forth in claim 1, and detecting the amplified product.
10. The method of 9, wherein the method further comprises detecting an amphotropic MLV virus using a primer pair and probe selected from:

SEQ ID NOs: 13, 14 and 15; and
SEQ ID NOs: 16, 17 and 18;
and/or the method further comprises detecting a cytosine deaminase gene using a primer and probe of SEQ ID NOs: 19, 20, and 21.

Patentansprüche

1. Zusammensetzung, die ein Primer-Paar und eine Sonde umfasst, ausgewählt aus der Gruppe bestehend aus:

(a) SEQ ID Nr. 7, 8 und 9; und
(b) SEQ ID Nr. 10, 11 und 12.
2. Verfahren zum Bestimmen des viralen Inhalts in einer *In-vitro*-Probe von einem Subjekt, das kurz vor einer retroviralen Genzuführungstherapie mittels eines MLV-bezogenen Virus steht oder sich in einer solchen befindet, das Folgendes beinhaltet:

Inkontaktbringen der Probe mit einem oder mehreren Primer-Paaren nach Anspruch 1 unter Bedingungen, die für eine Nukleinsäureamplifikation geeignet sind, um amplifizierte Produkte zu erhalten;
Inkontaktbringen der Probe mit einer oder mehreren Sonden, die mit dem amplifizierten Produkt hybridisieren;
Erkennen eines hybridisierten Produkts;
Anzeigen, dass das Subjekt viralen Inhalt hat, der ein MLV-bezogenes Virus umfasst.

3. Verfahren nach Anspruch 2, wobei das MLV-bezogene Virus ein in der Gentherapie benutzter rekombinanter retroviraler Vektor ist,
und/oder wobei das Verfahren vor einem Therapieplan ausgeführt wird, der einen Gentherapievektor umfasst,
und/oder wobei das Verfahren vor der Zuführung eines MLV-bezogenen retroviralen Vektors für Genzuführung ausgeführt wird,
und/oder wobei das Verfahren nach einer Zuführung eines MLV-bezogenen retroviralen Vektors für Genzuführung ausgeführt wird,
und/oder wobei das MLV-bezogene Virus ein 5' LTR, gag, pol, env Gene, eine regulatorische Domäne 3' zum env-Gen, verknüpft mit einem zuzuführenden heterologen Polynukleotid und ein 3' LTR und einen Promotor für die Expression in Säugetierzellen im 5'-LTR umfasst, wobei vorzugsweise die regulatorische Domäne eine interne Ribosomeneintrittsstelle (IRES) ist,
und/oder wobei vorzugsweise das heterologe Polynukleotid ein Polypeptid mit Cytosindeaminaseaktivität codiert,
und/oder wobei das Verfahren die Ausbreitung des MLV-bezogenen retroviralen Vektors überwacht, wobei das Verfahren vorzugsweise routinemäßig über den Verlauf von Jahren durchgeführt wird.

4. Verfahren zum Erkennen der Anwesenheit eines Virusmittels in einer Probe, das Folgendes beinhaltet:

Messen der Menge eines Polynukleotids in einer Probe mit einer quantitativen Polymerasekettenreaktion oder einem anderen Amplifikationsprozess, der ein Oligonukleotid-Primer-Paar und eine Sonde umfasst, ausgewählt aus:

- (i) SEQ ID Nr. 7, 8 und 9; und
- (ii) SEQ ID Nr. 10, 11 und 12.

5. Verfahren nach Anspruch 4, wobei das Polynukleotid DNA ist,
und/oder wobei das Polynukleotid RNA ist,
und/oder wobei die quantitative Polymerasekettenreaktion RT-qPCR ist,
und/oder wobei die Messung eine einzelne Kopie einer Virusmittel-bezogenen Nukleinsäure erkennt,
und/oder wobei das Virusmittel ein MLV-bezogenes Virus umfasst,
und/oder wobei die Probe Säugetiergewebe ist,
und/oder wobei die Probe Säugetierblut ist,
und/oder wobei das Virusmittel ein Gentherapievektor ist, wobei vorzugsweise der Gentherapievektor ein replikationskompetenter Vektor ist,
und/oder wobei das Verfahren vor einem Therapieplan ausgeführt wird, der eine Gentherapievektorbehandlung beinhaltet,
und/oder wobei das Verfahren zum Überwachen der Dosierung eines Therapieplans ausgeführt wird, der einen Gentherapievektor in einem Subjekt umfasst,
und/oder wobei der Gentherapievektor einen replikationskompetenten MLV-Vektor umfasst,
und/oder wobei das Verfahren nach einem Therapieplan ausgeführt wird, der einen Gentherapievektor umfasst,
und/oder wobei das Verfahren zum Überwachen der Dosierung eines Therapieplans durchgeführt wird, der einen Gentherapievektor umfasst,
und/oder wobei die Probe von einem menschlichen Subjekt ist.

6. Kit zum Ausführen des Verfahrens nach einem der Ansprüche 2 bis 5, wobei der Kit ein Primer-Paar gemäß Definition in Anspruch 1 umfasst.

7. Verfahren nach einem der Ansprüche 2 bis 5 zum Erkennen von <100 Kopien von MLV-bezogener DNA in einer Probe, extrahiert von festen histopathologischen Sektionen,
und/oder zum Erkennen von <100 Kopien von MLV-bezogener RNA in einer Probe, extrahiert von festen histopathologischen Sektionen.

8. Verfahren nach einem der Ansprüche 2 bis 5, wobei das Verfahren ein Verfahren zum Erkennen von MLV-bezogenen Viren in Plasma oder Serum von einem Säugetierwirt ist.

9. Verfahren zum Untersuchen eines Blutvorrats oder einer Gewebebank auf Infektion mit einem MLV-bezogenen Virus, umfassend eine Amplifikationsreaktion an dem Blutvorrat oder der Gewebebank unter Verwendung eines Primer-Paares nach Anspruch 1, und Erkennen des amplifizierten Produkts.

10. Verfahren nach Anspruch 9, wobei das Verfahren ferner das Erkennen eines amphotropen MLV-Virus mit einem Primer-Paar und einer Sonde beinhaltet, ausgewählt aus:

SEQ ID Nr. 13, 14 und 15; und

SEQ ID Nr. 16, 17 und 18;

und/oder das Verfahren ferner das Erkennen eines Cytosindeaminasegens mit einem Primer und einer Sonde von SEQ ID Nr. 19, 20 und 21 beinhaltet.

Revendications

1. Composition comprenant une paire d'amorces et une sonde sélectionnées dans le groupe consistant en les suivantes :

(a) SEQ ID NO : 7, 8 et 9 ; et

(b) SEQ ID NO : 10, 11 et 12.

2. Méthode de détermination *in vitro* du contenu viral dans un échantillon obtenu chez un sujet sur le point de recevoir ou en train de recevoir une thérapie par transfert d'un gène rétroviral qui utilise un virus apparenté au virus de la leucémie murine (MLV), comprenant :

la mise en contact de l'échantillon avec une ou plusieurs paires d'amorces telles que définies à la revendication 1 dans des conditions qui se prêtent à une amplification de l'acide nucléique pour obtenir des produits amplifiés ;

la mise en contact de l'échantillon avec une ou plusieurs sondes qui s'hybrident au produit amplifié ;

la détection d'un produit hybridé ;

indiquant que le sujet présente un contenu viral qui comprend un virus apparenté au MLV.

3. Méthode de la revendication 2, où le virus apparenté au MLV est un vecteur rétroviral recombinant utilisé en thérapie génique ;

et/ou où la méthode est effectuée avant un schéma de traitement qui comprend un vecteur de thérapie génique ;

et/ou où la méthode est effectuée avant l'administration d'un vecteur rétroviral apparenté au MLV servant au transfert d'un gène ;

et/ou où la méthode est effectuée après l'administration d'un vecteur rétroviral apparenté au MLV servant au transfert d'un gène ;

et/ou où le virus apparenté au MLV comprend une LTR 5', des gènes gag, pol et env, un domaine de régulation 3' du gène env lié à un polynucléotide hétérologue à transférer et une LTR 3' et, dans la LTR 5', un promoteur de l'expression dans des cellules mammifères, de préférence où le domaine de régulation est un site d'entrée interne des ribosomes (IRES) ;

et/ou de préférence où le polynucléotide hétérologue code pour un polypeptide ayant une activité cytosine désaminase ;

et/ou où la méthode surveille la propagation du vecteur rétroviral apparenté au MLV, de préférence où la méthode est effectuée en routine au cours des années de cure.

4. Méthode de détection de la présence d'un agent viral dans un échantillon comprenant :

la mesure de la quantité d'un polynucléotide dans un échantillon en utilisant une réaction en chaîne par polymérase quantitative ou un autre processus d'amplification comprenant une paire d'amorces oligonucléotidiques et une sonde sélectionnées parmi les suivantes :

(i) SEQ ID NO : 7, 8 et 9 ; et

(ii) SEQ ID NO : 10, 11 et 12.

5. Méthode de la revendication 4, où le polynucléotide est un ADN ;

et/ou où le polynucléotide est un ARN ;

et/ou où la réaction en chaîne par polymérase quantitative est une réaction en chaîne par polymérase quantitative en temps réel (RT-qPCR) ;

et/ou où la mesure détecte une copie unique d'un acide nucléique associé à un agent viral ;

et/ou où l'agent viral comprend un virus apparenté au MLV ;

et/ou où l'échantillon est un tissu mammifère ;

et/ou où l'échantillon est le sang d'un mammifère ;

et/ou où l'agent viral est un vecteur de thérapie génique, de préférence où le vecteur de thérapie génique est un vecteur compétent pour la réplication ;

et/ou où la méthode est effectuée avant un schéma de traitement qui comprend un traitement par un vecteur de thérapie génique ;

et/ou où la méthode est effectuée pour surveiller la posologie d'un schéma de traitement qui comprend un vecteur de thérapie génique chez un sujet ;

et/ou où le vecteur de thérapie génique comprend un vecteur MLV compétent pour la réplication ;

et/ou où la méthode est effectuée suite à un schéma de traitement qui comprend un vecteur de thérapie génique ;

et/ou où la méthode est effectuée pour surveiller la posologie d'un schéma de traitement qui comprend un vecteur de thérapie génique ;

et/ou où l'échantillon est obtenu chez un sujet humain.

6. Kit permettant de réaliser la méthode de l'une quelconque des revendications 2 à 5, le kit comprenant une paire d'amorces telles que définies à la revendication 1.

7. Méthode de l'une quelconque des revendications 2 à 5 pour la détection de <100 copies d'un ADN apparenté au MLV dans un échantillon extrait à partir de coupes histopathologiques fixées ;

et/ou pour la détection de <100 copies d'un ARN apparenté au MLV dans un échantillon extrait à partir de coupes histopathologiques fixées.

8. Méthode de l'une quelconque des revendications 2 à 5, où la méthode est une méthode de détection de virus apparentés au MLV dans le plasma ou le sérum d'un hôte mammifère.

9. Méthode de triage d'une banque de sang ou d'une banque de tissus pour le dépistage d'une infection par un virus apparenté au MLV, qui comprend une réaction d'amplification effectuée sur la banque de sang ou la banque de tissus en utilisant une paires d'amorces telles que définies à la revendication 1 et la détection du produit amplifié.

10. Méthode de la revendication 9, où la méthode comprend en outre la détection d'un virus de type MLV amphotrope en utilisant une paire d'amorces et une sonde sélectionnées parmi les suivantes :

SEQ ID NO: 13, 14 et 15; et

SEQ ID NO : 16, 17 et 18 ;

et/ou la méthode comprend en outre la détection d'un gène de la cytosine désaminase en utilisant une amorce et une sonde des SEQ ID NO : 19, 20 et 21.




FIGURE 1




FIGURE 2




FIGURE 3




FIGURE 4

Query	230	ACAGTTCCCGCCTCCGTCTGAATTTTGGCTTTCGG- <u>TTTGGGACCGAAGCCGCGCCGCGC</u>	288
Sbjct	486	ACACTTCCCGCCCCCGTCTGAATTTTGGCTTTCGGTTTTACG-CCGAAACCGCGCCGCGC	544
Query	289	GTCTTG-TCTG--CTGCAG--CATC-GTTCTGT-GTT-GTCTCTGTCTGACTGTGTTTCT	340
Sbjct	545	GTC-TGATTTGTTTTGTTGTTCTTCTGTCT-TCGTTAGT-T-T-TCT-TCTGTCTTTAA	598
Query	341	GTATTTGTCTGAGAAT-ATGGGCCAGACTGTTACCACTCCCT-TAAGTTTGACCTTAG--	396
Sbjct	599	GT-GTTCTC-GAG-ATCATGGGACAGACCGTAACCTAC-CCCTCTGAGTCTAACCTT-GCA	653
Query	397	GTCACCTGGAAAGATGTCGAGCGGATCGC-TCACAACCAGTCGGTAGATGTCAAGAAGAGA	455
Sbjct	654	G-CACTGGGGAGATGTCCAGCGCATTGCATC-CAACCAGTCTGTGGATGTCAAGAAGAGG	711
Query	456	CGTTGGGTACCTTCTGCTCTGCAGAATGGCCAAACCTTT-AACGTCGGATGGCCGCGA-G	513
Sbjct	712	CGCTGGGTACCTTCTGTTCCGCCGAATGGCCAA-CTTTCAATGTAGGATGGCCTC-AGG	769
Query	514	ACGGCACCTTTAACCAGACCTCATCACCAGGTTAAGATCAAG-GTCTTTTCACCTGG-	571
Sbjct	770	ATGGTACTTTTAATTTAGGTGTTATCTCTCAGGTCAAG-TCTAGAGTGTGTTTGTCTGGT	828
Query	572	CCCGCATGGACACCCAGACCAGGTCCCCTACATCGTGACCTGGGAAGC-CTTGGCTTTTG	630
Sbjct	829	CCC-CACGGACACCCGGATCAGGTCCCATATATCGTCACCTGGGAGGCACTT-GCCTATG	886
Query	631	ACCCCCCTCCCTGGGTCAAGCCCTTTGTACAC-CCTAAGCCTCCGCCTCCTCTTCCTCCA	689
Sbjct	887	ACCCCCCTCCGTGGGTCAAACCGTTTGT-CTCTCCTAA-AC-CC-CCTCCT-TTAC-CGA	940
Query	690	TCCGC-CCCGTCTCTCCC-CCTTGAACCTCCTCGTTCGACC-CCG-CCTCGATCCT-CCC	744
Sbjct	941	-CAGCTCCCGTC-CTCCCGCC-CGGTCCTTCT-GCGCAACCTCCGTCC-CGAT-CTGCCC	994
Query	745	TTTATCCAGCCCTCACTCCTTCTCT-AGGCGCCAAA---CCTAAACCTCAAGTTCTTTCT	800
Sbjct	995	TTTACCCTGCCCTTACCCCTCTATAAAG-TCCAAACCTCCTAAGCCCCAGGTTCTCCCT	1053
Query	801	GACAGTGGGGGGCCGCTCATCGACCTACTTACAGAAGACCCCCCGCCTTATAGGGACCCA	860
Sbjct	1054	GATAGCGGCGGACCTCTCATTGACCTTCTCACAGAGGATCCCCCGCC--GTACGGA-GCA	1110
Query	861	-AGACC-ACCCCTTCCGACAGGGACGGA-AATGGTGGAGAAGC-G--ACC-CCTGCGGG	913
Sbjct	1111	CA-ACCTTCCTCC-TCTGCCAGGGA-GAACAATGAAGAAGAGGCGGCCACCACCT-C-CG	1165
Query	914	AGAGGCACCGGACCCC-TC-CCCAATGGCATCTCGCCTACGTGGGAGACGG-GAG-CCC-	968
Sbjct	1166	AG-GTTTCC--CCCCCTTCTCCC-ATGGTGTCTCGACTGCG-GGGA-A-GGAGAGACCCT	1218
Query	969	CCTGTGGCCGACTCCACTACCTCGCAGGCATTCCCCCTCCGCGCAGGAGGAAACGGACAG	1028
Sbjct	1219	CCCGCAGCGGACTCCACCACCTCCCAGGCATTCCCACTCCGCATGGGGGGAGATGGCCAG	1278

FIGURE 5

Query	1029	CTTCAATACTGGCCGTTCTCCTCTTCTGACCTT-TACAAC <u>TGGAAAAATAATAACCCCTTC</u>	1087
Sbjct	1279	CTTCAGTACTGGCCGTTTTCCTCCTCTGA-TTTATATAATTGGAAAAATAATAACCCCTTC	1337
Query	1088	<u>TTTTCTGAAGATCCAGGTAAACTGACAG</u> -CTCTGATCGAGTCTGTTCTCATCACCCATC	1146
Sbjct	1338	CTTTTCTGAAGATCCAGGTAAATTGACGGCCT-TGATTGAGTCCGTCTCATCACCCACC	1396
Query	1147	<u>AGCCACCTGGGACGACTGTCAGCAGCTGTTGGGGACTCTGCTGACCGGAGAAGAAAAAC</u>	1206
Sbjct	1397	AGCCACCTGGGACGACTGTCAGCAGTTGTTGGGGACCCTGCTGACCGGAGAAGAAAAAGC	1456
Query	1207	AACGGGTGCTCTTAGAGGCTAGAAAGGCGGTGCGGGGCGATGATGGGCGCCCCACTCAAC	1266
Sbjct	1457	AGCGGGTGCTCCTAGAGGCTAGAAAGGCAGTCCGGGGCAATGATGGACGCCCACTCAGT	1516
Query	1267	TGCCCCAATGAAGTCGATGCCGCTTTTCCCTCGAGCGCCAGACTGGGATTACACCAC-C	1325
Sbjct	1517	TGCCTAATGAAGTCAATGCTGCTTTTCCCTTGAGCGCCCCGATTGGGATTACACCACTA	1576
Query	1326	<u>CAGGCAGGTAGGAACCACCTAGTCCACTATCGCCAGTTGCTCCTAGCGGGTCTCCAAAAAC</u>	1385
Sbjct	1577	CA-GAAGGTAGGAACCACCTAGTCCTCTACCGCCAGTTGCTCTTAGCGGGTCTCCAAAAAC	1635
Query	1386	<u>GCGGGCAGAAGCCCCACCAATTTGGCCAAGGTAAAAGGAATAACACAAGGG-CCCAATGA</u>	1444
Sbjct	1636	GCGGGCAGGAGCCCCACCAATTTGGCCAAGGTAAAAGGGATAAC-CCAGGGACCTAATGA	1694
Query	1445	GTCTCCCTCGGCCTTCCTAGAGAGACTTAAGGAAGCCTATCGCAGGTACACTCCTTATGA	1504
Sbjct	1695	GTCTCCCTCAGCCTTTTATAGAGAGACTCAAGGAGGCCTATCGCAGGTACACTCCTTATGA	1754
Query	1505	<u>CCCTGAGGACCCAGGGCAAGAACTAATGTGTCTATGTCTTTTCAATTTGGCAGTCTGCCCC</u>	1564
Sbjct	1755	CCCTGAGGACCCAGGGCAAGAAACCAATGTGTCCATGTCATTTCATCTGGCAGTCTGCCCC	1814
Query	1565	AGACATTGGGA-GAAAGTTAGAGAGGTTAGAAGATTTAAAAAACAAGACGCTT-GGAGAT	1622
Sbjct	1815	GGATATCGG-ACGAAAGTTAGAGCGGTTAGAAGATTTAAAGAGCAAGAC-CTTAGGAGAC	1872
Query	1623	TTGGTTAGAGAGGCAGAAAAGATCTTTAATaaacgaagaaacccctggaaggaatgagatgaa	1682
Sbjct	1873	TTAGTGAGGGAAGCTGAAAAGATCTTTAATAAGCGAGAAACCCCGGAAGAAAGAGAGGAA	1932
Query	1683	<u>cgtatcagtgagagaaaacagaggaanagaggaacgcctatggaacagaggaatgagcagaa</u>	1742
Sbjct	1933	CGTATCAGGAGAGAAATAGAGGAAAAAGAAGAACGCCGTAGGGCAGAGGATGAGCAGAGA	1992
Query	1743	<u>gagaaagaaaggaatTCGTAGGAGACATAGAGAGATGAGCAAGCTATTGGCCACTGTCGTT</u>	1802
Sbjct	1993	GAGAGAGAAAGGGACCGCAGAAGACATAGAGAGATGAGCAAGCTCTTGGCCACTGTAGTT	2052
Query	1803	AGTGGACAGAAACAGGATAGACAGGGAGGAGAACGAAGGAGGTCCCAACTCGATC-GCGA	1861
Sbjct	2053	ATTGGTCAGAGACAGGATAGACAGGGGGGAGAGCGGAGGAGGCCCAACTTGATAAG-GA	2111

FIGURE 5 (cont'd)

EP 2 913 402 B1

Query	1862	CCAGTGTGCTCTACTGCAAAGAAAAGGGGCACTGGGCTAAAGATTGTCCCAA-GAAACCAC	1920
Sbjct	2112	CCAATGCGCTACTGCAAAGAAAAGGGACACTGGGCTAAGGACTG-CCCAAAGAAGCCAC	2170
Query	1921	GAGGACCTCGGGGACCAAGACCCCAGACCTCCCTCCTGACCCTAGATGACTAGGGAGGTC	1980
Sbjct	2171	GAGGGCCCCGAGGACCGAGGCCCCAGACCTCCCTCCTGACCTTAGGTGACTAGGGAGGTC	2230
Query	1981	AGGGTCAGGAGGCGGCTGTAACCCAGGATAACCCCTCAAAGTCGGGGGGCAACCCGTCA	2040
Sbjct	2231	AGGGTCAGGAGCCCCCCCCTGAACCCAGGATAACCCCTCAAAGTCGGGGGGCAACCCGTCA	2290
Query	2041	CCTTCCTGGTAGATACTGGGGCCCCAACACTCCGTGCTGACCCAAAATCCTGGACCCCTAA	2100
Sbjct	2291	CCTTCCTGGTAGATACTGGGGCCCCAACACTCCGTGCTGACCCAAAATCCTGGACCCCTAA	2350
Query	2101	GTGATAAGTCTGCCTGGGTCCAAGGGGCTACTGGAGGAAAGCGGTATCGCTGGACCACGG	2160
Sbjct	2351	GTGACAAGTCTGCCTGGGTCCAAGGGGCTACTGGAGGAAAGCGGTATCGCTGGACCACGG	2410
Query	2161	ATCGCAAAGTACATCTAGCTACCGGTAAGGTCACCCACTCTTTCCTCCATGTACCAGACT	2220
Sbjct	2411	ATCGCAAAGTACATCTGGCTACCGGTAAGGTCACCCACTCTTTCCTCCATGTACCAGACT	2470
Query	2221	GTCCCTATCCTCTGTTAGGAAGAGATTGCTGACTAAACTAAAAGCCCAAATCCACTTTG	2280
Sbjct	2471	GCCCTATCCTCTGCTAGGAAGAGACTTGCTGACTAAACTAAAAGCCCAAATCCACTTTG	2530
Query	2281	AGGGATCAGGAGCTCAGGTTATGGGACCAATGGGGCAGCCCCTGCAAGTGTGACCCTAA	2340
Sbjct	2531	AGGGATCAGGAGCTCAGGTTGTGGGACCGATGGGACAGCCCCTGCAAGTGTGACAGTAA	2590
Query	2341	ATATAGAAGATGAGCATCGGCTACATGAGACCTCAAAAGAGCCAGATGTTTCTCTAGGGT	2400
Sbjct	2591	ACATAGAAGATGAGTATTGGCTACATGATACCAGGAAAGAGCCAGATGTTCTCTAGGGT	2650
Query	2401	CCACATGGCTGTCTGATTTTCCT-CAGGCCTGGGCGGAAACCGGGGGCATGGGACTGGCA	2459
Sbjct	2651	CCACATGGCTTTCTGATTT-CCTTCAGGCCTGGGCGGAAACCGGGGGCATGGGACTGGCA	2709
Query	2460	GTTCGCCAAGCTCCTCTGATCATACCTCTGAAAGCAACCTCTACCCCCGTGTCCATAAAA	2519
Sbjct	2710	GTTCGCCAAGCTCCTCTGATCATACCTCTGAAGGCAACCTCTACCCCCGTGTCCATAAAA	2769
Query	2520	CAATACCCCATGTACAGAAGCCAGACTGGGGATCAAGCCCCACATACAGAGACTGTTG	2579
Sbjct	2770	CAATACCCCATGTACAGAAGCCAGACTGGGGATCAAGCCCCACATACAGAGGCTGTTG	2829
Query	2580	GACCAGGGAATACTGGTACCCTGCCAGTCCCCCTGGAACACGCCCCTGCTACCCGTTAAG	2639
Sbjct	2830	GACCAGGGAATACTGGTACCCTGCCAGTCCCCCTGGAACACGCCCCTGCTACCCGTTAAG	2889
Query	2640	AAACCAGGGACTAATGATTATAGGCCTGTCCAGGATCTGAGAGAAGTCAACAAGCGGGTG	2699
Sbjct	2890	AAACCAGGGACTAATGATTATAGGCCTGTCCAGGATCTGAGAGAAGTCAACAAGCGGGTG	2949

FIGURE 5 (cont'd)

Query	2700	<u>GAAGACATCCACCCCACCGTGCCCAACCCTTACAACCTCTTGAGCGGGCTCCCACCGTCC</u>	2759
Sbjct	2950	GAAGACATCCACCCCACCGTGCCCAACCCTTACAACCTCTTGAGCGGGCTCCCACCGTCC	3009
Query	2760	<u>CACCAGTGGTACACTGTGCTTGATTAAAGGATGCCTTTTTCTGCCTGAGACTCCACCCC</u>	2819
Sbjct	3010	CACCAGTGGTACACTGTGCTTGATTAAAGGATGCCTTTTTCTGCCTGAGACTCCACCCC	3069
Query	2820	<u>ACCAGTCAGCCTCTCTTCGCCTTTGAGTGGAGAGATCCAGAGATGGGAATCTCAGGACAA</u>	2879
Sbjct	3070	ACCAGTCAGCCTCTCTTCGCCTTTGAGTGGAGAGATCCAGAGATGGGAATCTCAGGACAA	3129
Query	2880	<u>TTGACCTGGACCAGACTCCCACAGGGTTTCAAAAACAGTCCCACCCTGTTTGATGAGGCA</u>	2939
Sbjct	3130	CTGACCTGGACCAGACTCCCACAGGGTTTCAAAAACAGTCCCACCCTGTTTGATGAGGCA	3189
Query	2940	<u>CTGCACAGAGACCTAGCAGACTTCCGGATCCAGCACCCAGACTTGATCCTGCTACAGTAC</u>	2999
Sbjct	3190	CTGCACAGAGACCTAGCAGATTTCGGATCCAGCACCCAGACTTGATCCTGCTACAGTAC	3249
Query	3000	<u>GTGGATGACTTACTGCTGGCCGCCACTTCTGAGCTAGACTGCCAACAAGGTACTCGGGCC</u>	3059
Sbjct	3250	GTGGATGACTTACTGCTGGCCGCCACTTCTGAGCAAGACTGCCAACGAGGTACTCGGGCC	3309
Query	3060	<u>CTGTTACAAACCCTAGGGAACCTCGGGTATCGGGCCTCGGCCAAGAAAGCCCAAATTTGC</u>	3119
Sbjct	3310	CTATTACAAACCCTAGGGAACCTCGGGTATCGGGCCTCGGCCAAGAAAGCCCAAATTTGC	3369
Query	3120	<u>CAGAAACAGGTCAAGTATCTGGGGTATCTTCTAAAAGAGGGTCAGAGATGGCTGACTGAG</u>	3179
Sbjct	3370	CAGAAACAGGTCAAGTATCTGGGGTATCTCCTAAAAGAGGGACAGAGATGGCTGACTGAG	3429
Query	3180	<u>GCCAGAAAAGAGACTGTGATGGGGCAGCCACTCCGAAGACCCCTCGACAACCTAAGGGAG</u>	3239
Sbjct	3430	GCCAGAAAAGAGACTGTGATGGGGCAGCCACTCCGAAGACCCCTCGACAACCTAAGGGAG	3489
Query	3240	<u>TTCTTAGGGACGGCAGGCTTCTGTGCCTCTGGATCCCTGGGTTTGCAGAAATGGCAGCC</u>	3299
Sbjct	3490	TTCTTAGGGACGGCAGGCTTCTGTGCCTCTGGATCCCTGGGTTTGCAGAAATGGCAGCC	3549
Query	3300	<u>CCCTTGTACCCTCTACCAAAACGGGGACTCTGTTTAATTGGGGCCCAGACCAACAAAAG</u>	3359
Sbjct	3550	CCCTTGTACCCTCTTACCAAAACGGGGACTCTGTTTAATTGGGGCCCAGACCAGCAAAAG	3609
Query	3360	<u>GCCTATCAAGAAATCAAGCAAGCTCTTCTAACTGCCCCAGCCCTGGGGTTGCCAGATTTG</u>	3419
Sbjct	3610	GCCTATCAAGAAATCAACAGGCTCTTCTAACTGCCCCGCGCTGGGATTGCCAGATTTG	3669
Query	3420	<u>ACTAAGCCCTTTGAACTCTTTGTGACGAGAAGCAGGGCTACGCCAAAGGTGTCCTAACG</u>	3479
Sbjct	3670	ACTAAGCCCTTTGAACTCTTTGTGACGAGAAGCAGGGCTACGCCAAAGGCGTCCTAACG	3729
Query	3480	<u>CAAAAACCTGGGACCTTGGCGTCGGCCGGTGGCCTACCTGTCCAAAAAGCTAGACCCAGTA</u>	3539
Sbjct	3730	CAAAAACCTGGGACCTTGGCGTCGGCCGTGGCCTACCTGTCCAAAAAGCTAGACCCAGTG	3789

FIGURE 5 (cont'd)

Query	3540	<u>GCAGCTGGGTGGCCCCCTTGCCTACGGATGGTAGCAGCCATTGCCGTACTGACAAAGGAT</u>	3599
Sbjct	3790	GCAGCTGGGTGGCCCCCTTGCCTACGGATGGTAGCAGCCATTGCCGTCTGACAAAAAAT	3849
Query	3600	GCAGGCAAGCTAACCATGGGACAGCCACTAGTCATTCTGGCCCCCATGCAGTAGAGGCA	3659
Sbjct	3850	GCAGGCAAGCTAACTATGGGACAGCCGCTAGTCATTCTGGCCCCCATGCGGTAGAAGCA	3909
Query	3660	CTAGTCAAACAACCCCCGACCGTGGCTTTCCAACGCCCGGATGACTCACTATCAGGCC	3719
Sbjct	3910	CTGGTCAAACAACCCCCGTGACCGTGGCTATCCAATGCCCGCATGACCCACTATCAGGCA	3969
Query	3720	TTGCTTTTGGACACGGACCGGGTCCAGTTCGGACCGGTGGTAGCCCTGAACCCGGCTACG	3779
Sbjct	3970	ATGCTCCTGGATACAGACCGGGTTCAGTTCGGACCGGTGGTGGCCCTCAACCCGGCCACC	4029
Query	3780	CTGCTCCCACTGCCTGAGGAAGGGCTGCAACACA--ACTGCCTTGATATCCTGGCCGAAG	3837
Sbjct	4030	CTGCTCCCCCTACCGGAA-AAGGAA-GCCCCCATGACTGCCTCGAGATCTTGGCTGA-G	4086
Query	3838	CC-CACGGAACCCGACCCGACCTAACGGACCAGCCGC-TCCCAGACGCCGACCACACCTG	3895
Sbjct	4087	ACGCACGGAACCAGACCGGACCTCACGGACCAGCC-CATCCCAGACGCTGATTACACTTG	4145
Query	3896	GTACACGGATGGAAGCAGTCT-CTTACAAGAGGGACAGCGTA-AGGCGGGAGCTGCGGTG	3953
Sbjct	4146	GTACACAGATGGAAGCAG-CTTCCTACAAGAAGGACAACGGAGAG-CTGGAGCAGCGGTG	4203
Query	3954	ACCACCGAGACCGAGGTAATCTGGGCTAAAGCCCTGCCAGCCGGGACATCCGCTCAGCGG	4013
Sbjct	4204	ACTACTGAGACCGAGGTAATCTGGGCGAGGGCTCTGCCGGCTGGAACATCCGCCCAACGA	4263
Query	4014	<u>GCTGAACTGATAGCACTCACCCAGGCCCTAAAGATGGCAGAAGGTAAGAAGCTAAATGTT</u>	4073
Sbjct	4264	GCCGAACTGATAGCACTCACCCAAGCCTTAAAGATGGCAGAAGGTAAGAAGCTAAATGTT	4323
Query	4074	TATACTGATAGCCGTTATGCTTTTGGCTACTGCCCATATCCATGGAGAAATATACAGAAGG	4133
Sbjct	4324	TACACTGATAGCCGCTATGCCTTCGCCACGGCCCATGTCCATGGAGAAATATATAGGAGG	4383
Query	4134	CGTGGGTTGCTCACATCAGAAGGCAAAGAGATCAAAAATAAGA-CGAGATCTTGGCCCT	4192
Sbjct	4384	CGAGGGTTGCTGACCTCAGAAGGCAGAGAAATTAATAACAA-GAACGAGATCTTGGCCTT	4442
Query	4193	ACTAAAAGCCCTCTTTCTGCCCCAAAAGACTTAGCATAATCCATTGTCCAGGACATCAAAA	4252
Sbjct	4443	GCTAAAAGCTCTCTTTCTGCCCCAACGACTTAGTATAATCACTGTCCAGGACATCAAAA	4502
Query	4253	GGGACACAGCGCCGAGGCTAGAGGCAACCGGATGGCTGACCAAGCGGCCCGAAAGGCAGC	4312
Sbjct	4503	AGGAAACAGTGCTGAGGCCAGAGGCAACCGTATGGCAGATCAAGCAGCCCGAGAGGCAGC	4562
Query	4313	CATCACAGAG-ACTCC-AGACACCTCTACCCTCCTCATAGAAAATTCATCACCC-TACAC	4369
Sbjct	4563	CATGA-AG-GCAGTTCTAGAAAACCTCTACACTCCTCATAGAGGACTCAAC-CCCGTATAC	4619

FIGURE 5 (cont'd)

EP 2 913 402 B1

Query	4370	-C-TCAGAACATTTTCATTACACAGTGACTGATAT-AAAGGACCTAACC--AAGTTGGGG	4424
Sbjct	4620	GCCTCCC--CATTTCCATTACACCGAAACAGATCTCAAAGAC-TA-CGGGAACT-GGGA	4674
Query	4425	GCCATTTATGATAAAACAAA-GAAGTATTGGGTC-TACCAAGGAAAACCTGTGATGCCTG	4482
Sbjct	4675	GCCACCTACAATCAGACAAAAGGA-TATTGGGTCTTAC-AAGGCAAACCTGTGATGCCCCG	4732
Query	4483	ACCAGTTTACT-TTTGAATTATTAGACTTCTTCATCAG-CTGACTCACCTCAGCTTCTC	4540
Sbjct	4733	ATCAGTCCG-TGTTTGAAGTGTAGACTCCCTACA-CAGACTCACCCATCTGAGCC-CTC	4789
Query	4541	AAAA-ATGAAGGCTCTCCTAGA--GAGAAGCCACAGTCCCTACTACATGCTGAACCGGGA	4597
Sbjct	4790	AAAAGATGAAGGCACTCCTCGACAGAGAAGA-A-AGCCCCTACTACATGTTAAACCGGGA	4847
Query	4598	TC-GAAC-ACTCAAAAATATCACTGAGACCTGCAAAGCTTGTGCACAAGTCAACGCCAGC	4655
Sbjct	4848	-CAGAATA-TCCAGTATGTGACTGAGACCTGCACCGCCTGTGCCAAGTAAATGCCAGC	4905
Query	4656	AA-GTCTGCCGTTAAA-CAGGGAA-CTAGGGTCCGCGGGCATCGGCCCGGCACTCATTGG	4712
Sbjct	4906	AAAGCCAAAA-TTGGGGCAGGGGTGCGAG--TACGCGGACATCGGCCAGGCACCCATTGG	4962
Query	4713	GA-GATCGATTTACCGA-GATAAAGCCCGGATTGTATGGCTATAAATATCTTCTAGTTT	4770
Sbjct	4963	GAAGTT-GATTTACGGAAG-TAAAGCCAGGACTGTATGGGTACAAGTACCTCCTAGTGT	5020
Query	4771	TTATAGATACCTTTTCTGGCTGGATAGAAGCCTTCCCAACCAAGAAAGAAACCGCCAAGG	4830
Sbjct	5021	TTGTAGACACCTTCTCTGGCTGGGTAGAGGCATTCCCGACCAAGCGGGAAACTGCCAAGG	5080
Query	4831	TCGTAACCAAGAAGCTACTAGAGGAGATCTTCCCCAGGTTCCGGCATGCCTCAGGTATTGG	4890
Sbjct	5081	TCGTGTCCAAAAAGCTGTTAGAAGACATTTTTCCGAGATTGGAAATGCCGCAGGTATTGG	5140
Query	4891	GAACTGACAATGGGCCTGCCTTCGTCTCCAAGGTGAGTCAGACAGTGGCCGATCTGTTGG	4950
Sbjct	5141	GATCTGATAACGGGCCTGCCTTCGCCTCCAGGTAAGTCAGTCAGTGGCCGATTTACTGG	5200
Query	4951	GGATTGATTGGAAATTACATTGTGCATACAGACCCCAAAGCTCAGGCCAGGTAGAAAGAA	5010
Sbjct	5201	GGATCGATTGGAAGTTACATTGTGCTTATAGACCCAGAGTTCAGGACAGGTAGAAAGAA	5260
Query	5011	<u>TGAATAGAACCATCAAGGAGACTTTAACTAAATTAACGCTTGCAACTGGCTCTAGAGACT</u>	5070
Sbjct	5261	TGAATAGACAATTAAGGAGACTTTGACCAAATTAACGCTTGCACTGGCACTAGAGACT	5320
Query	5071	<u>GGGTGCTCCTACTCCCTTAGCCCTGTACCGAGCCCGCAACACGCCGGGCCCCCATGGCC</u>	5130
Sbjct	5321	GGGTACTCCTACTCCCTTAGCCCTCTACCGAGCCCGGAATACTCCGGGCCCCACGGAC	5380
Query	5131	<u>TCACCCCATATGAGAT-CTTATATGGGGCACCCCGCCCTTGTA</u> AACTTCCCTGACCCT	5189
Sbjct	5381	TGACTCCGTATGAAATTCTG-TATGGGGCACCCCGCCCTTGTCCAATTTTCATGATCCT	5439

FIGURE 5 (cont'd)

Query	5190	GACATGACAAGAGTTA-CTAACAGCC CCTCTCTCCAAGCTCACTTACAGGGCTCTCTACTT	5248
Sbjct	5440	GAAATGTCAA-AGTTAACTAATAGTCCCTCTCTCCAAGCTCACTTACAGGCCCTCCAAGC	5498
Query	5249	AGTCCAGCACGAAGTCTGGA-GACCTCTGGCGGCAGCCTACCAAGAACAACCTGGACC-GA	5306
Sbjct	5499	AGTACAACAAGAGGTCTGGAAG-CCGCTGGCCGCTGCTTATCAGGACCAGCTAGATCAG-	5556
Query	5307	CCGGTGGTACCTCACCC-TTACCGAGTCGGCGACACAGTGTGGGTCCGCCGACACCAGAC	5365
Sbjct	5557	CCAGTGATACCACACCCCTT-CCGTGTCGGTGACGCCGTGTGGGTACGCCGGCACCAGAC	5615
Query	5366	TAAGAACCTAGA ACCTCGCTGGAAAGGACCTTACACAGTCTGTGCTGACCACCCCCACCGC	5425
Sbjct	5616	TAAGAACTTAGAACCTCGCTGGAAAGGACCTTACACCGTCTGTGCTGACAACCCCCACCGC	5675
Query	5426	CCTCAAAGTAGACGGCATCGCAGCTTGGATACACGCCGCCACGTGAAGGCTGCCGAC	5483
Sbjct	5676	TCTCAAAGTAGACGGCATCTCTGCGTGGATACACGCCGCTCACGTAAAGGCGGC-GAC	5732
.			
.			
.			
.			
.			
Query	6436	GGGAACGGAAAAATAGGCTGCTAAACTTAGTAGACGGAGCCTACCAAGCCCTCAACCTCAC	6495
Sbjct	6603	GGG-ACGGGAGACAGGCTGCTAAACCTGGTAGAAGGAGCCTACCTAGCCCTCAACCTCAC	6661
Query	6496	<u>CAGTCCTGACAAAACCCAAGAGTGCTGGTTGTGTCTAGTAGCGGGACCCCCCTACTACGA</u>	6555
Sbjct	6662	CAGTCCCACAAAACCCAAGAGTGCTGGCTGTGTCTAGTATCGGGACCCCCCTACTACGA	6721
Query	6556	<u>AGGGGTTGCCGTCCTGGGTACCTACTCCAACCATACTCTGCTCCAGCCAACTGCTCCGT</u>	6615
Sbjct	6722	AGGGGTGGCCGTCCTAGGTACTTACTCCAACCATACTCTGCCCGGCTAACTGCTCCGT	6781
Query	6616	GGCCTCCCAACACAAGTTGACCCTGTCCGAAGTGACCGGACAGGGACTCTGCATAGGAGC	6675
Sbjct	6782	GACCTCCCAACACAAGCTGACCCTGTCCGAAGTGACCGGGCAGGGACTCTGCATAGGAGC	6841
Query	6676	<u>AGTTCCCAAAACACATCAGGCCCTATGTAATACCACCCAGACA-A-GCAGTCGA-GGGTC</u>	6732
Sbjct	6842	AGTTCCCAAAACCCATCAGGCCCTGTGTAATACCACCCAGA-AGACG-AG-CGACGGGTC	6898
Query	6733	CTATTATCTAGT-TGC-CCCTACAGGTACCATGTGGGCTTGTAGTACCGGGCTTACTCCA	6790
Sbjct	6899	CTACTATTTGGCCT-CTCCCGCCGGG-ACCATTTGGGCTTG CAGCACCGGGCTCACTCCC	6956
Query	6791	TG-C-ATCTCCACCACCATACTGAACCTTA-CCACTGATTATTGTGTTCTTGTGCAACTC	6847
Sbjct	6957	TGTCTATCT--ACTACTGTGCTTAAC-TTAACCACTGATTACTGTGTCCTGGTTGAACTC	7013

FIGURE 5 (cont'd)

Query	6848	TGGCCAAGAG-TCACCTATCATTCCCCAGCTATGTTTACGGCCTGTTTGAGA-GATCCA	6905
Sbjct	7014	TGGCCAA-AGGTAACCTACCACTCCCCTAATTATGTTTATGGCCAGTTTGAAAAGAA--A	7070
Query	6906	ACCGAC-ACAAAAGAGAACCGGTGTCGTTAACCCCTGGCCCTATTATTGGGTGGACTAACC	6964
Sbjct	7071	ACTAAATATAAAAAGAGAGCCGGTGTCATTAACCTCTGGCCCTGCTGTTGGGAGGACTTACT	7130
Query	6965	ATGGGGGGAATTGCCGCTGGAATAGGAACAGGGACTACTGCTCTAATGGCCACTCAG-CA	7023
Sbjct	7131	ATGGGCGGCATAGCTGCAGGAGTTGGAACAGGGACTACAGCCCTAGTGGCCAC-CAAACA	7189
Query	7024	ATTCCAGCAGCTCCAAGCCGC-AGTACAGG-ATGATCTCAGGGAGGTT-GAAAAATCAAT	7080
Sbjct	7190	ATTCGAGCAGCTCCAGGCAGCCA-TACATACA-GACCTTGGGGCC-TTAGAAAAATCAGT	7246
Query	7081	CTCTAACCTAGAAAAGTCTCTCACTTCCCTGTCTGAAGTTGTCCTACAGAATCGAAGGGG	7140
Sbjct	7247	CAGTGCCCTAGAAAAGTCTCTGACCTCGTTGTCTGAGGTGGTCTACAGAACCGGAGGGG	7306
Query	7141	CCTAGA-CTTGTTATTTCTAAAAGAAGGAGGGCTGTGTGCTGCTCTAAAAGAAGAATGTT	7199
Sbjct	7307	ATTAGATCTACTG-TTCCTAAAAGAAGGAGGATTATGTGCTGCCCTAAAAGAAGAATGCT	7365
Query	7200	GCTTCTATGCGGACCACACAGGAC-TAGTGAGAGACAGCATGGCCAAAT-TGAGAGAGAG	7257
Sbjct	7366	GTTTTTACGCGGACCACACTGG-CGTAGTAAGAGATAGCATGGC-AAAGCTAAGAGAAAG	7423
Query	7258	GCTTAATC-AGAGACAGAAACTGTTTGAGTCAACTCAAGGATGGTTTGAGGGACTGTTTA	7316
Sbjct	7424	G-TTAAACCAGAGACAAAATTGTTTCAATCAGGACAAGGGTGTTTGAGGGACTGTTTA	7482
Query	7317	<u>ACAGATCCCCTTGGTTTACCACCTTGATATCTACCATTATGGGACCC-CTCATTGTACTC</u>	7375
Sbjct	7483	ACAGGTCCCCTGGTTTACGACCCTGATATCCACCATTATGGG-CCCTCTGATAGTACTT	7541
Query	7376	CTAATGATTTTGCTCTTCGACCCTGCATTCTTAATCGATTAGTCCAATTTGTTAAAGAC	7435
Sbjct	7542	TTATTAATCCTACTCTTCGACCCTGTATTCTCAACCGCTTGGTCCAGTTTGTTAAAGAC	7601
Query	7436	AGGATATCAGTGGTCCAGGCTCTAGTTTTGACTCAACAATATCACCAGCTGAAGCCTATA	7495
Sbjct	7602	AGAATTTTCGGTAGTGCAGGCCCTGGTTCTGACCCAACAGTATCACCAACTCAAATCAATA	7661
Query	7496	GAGTAC-GA-GCCA-TAGA-TAA----AATAAAAGATTTTATTTAGTCTCCAGAAAAAGG	7547
Sbjct	7662	GA-TCCAGAAGA-AGTGGAATCACGTGAATAAAAGATTTTATTCAGTTTCCAGAAAGAGG	7719
Query	7548	<u>GGGGAATGAAAGACCCACCTGTA-GGTTTGGCAAGCTAGCTTA-AGTAACGCCATTTTG</u>	7605
Sbjct	7720	GGGGAATGAAAGACCCACCA-TAAGGCTTAGCACGCTAGCT-ACAGTAACGCCATTTTG	7777
Query	7606	<u>CAAGGCATGGAAAAATAC-ATAACTGAGAA-T-AGAGAAGTT-CA-G-A--TC-AAGGTC</u>	7656
Sbjct	7778	CAAGGCATGGAAAAGTACCAGAGCTGAGTTCTCA-A-AAGTTACAAGGAAGTTTAAAT-TA	7834

FIGURE 5 (cont'd)

EP 2 913 402 B1

[illegible]

FIGURE 5 (cont'd)




FIGURE 6




FIGURE 7

XMRV-gag
Standards

1-Stage qPCR_pUC57 XMRV gag Plasmid Standards			
pUC57 XMRV gag Input Copies/Reaction	N Rows	Mean Ct	Standard Deviation
1E0	2	34.88	0.23
1E1	3	30.37	0.16
1E2	3	27.07	0.03
1E3	3	23.27	0.14
1E4	3	20.39	0.09
1E5	3	16.63	0.03
1E6	3	14.43	0.07
1E7	3	12.51	0.04

FIGURE 8A

Controls

1-Stage qPCR_XMRV gag Controls				
Controls	N Rows	Mean Ct	Standard Deviation	Mean Copies/Reaction
22Rv1 Pos	3	14.63	0.03	870088.25
h.blood Neg	3	ND	ND	ND
NTC	3	ND	ND	ND

FIGURE 8B

Experimental (in human blood gDNA)

1-Stage qPCR_pUC57-XMRV gag Spiked Human Blood gDNA						
Sample (human blood gDNA)	pUC57-XMRV gag	N Rows	Mean Ct	Standard Deviation	Copies/ reaction	%Recovery of Copies/reaction
Neat	1E0	3	34.31	0.42	0.58	58.00
	1E1	3	31.04	0.15	6.09	60.90
	1E2	3	27.56	0.14	75.23	75.23
	1E3	3	23.83	0.11	1120.42	112.04
	1E4	3	21.46	0.1	6206.85	62.07
	1E5	3	17.48	0.19	110994.26	110.99
	1E6	3	15.23	0.08	564101.99	56.41
	1E7	3	12.91	0.07	3026602.41	30.27

FIGURE 8C

XMRV-env

Standards

1-Stage qPCR_pET28b XMRV env Plasmid Standards			
pET28b XMRV env Input Copies/Reaction	N Rows	Mean Ct	Standard Deviation
1E0	3	33.68	0.98
1E1	3	29.56	0.13
1E2	3	26.5	0.08
1E3	3	23.31	0.09
1E4	3	20.1	0.12
1E5	3	16.42	0.06
1E6	3	14.41	0.05
1E7	3	12.07	0.02

FIGURE 9A

Controls

1-Stage qPCR_XMRV env Controls				
Controls	N Rows	Mean Ct	Standard Deviation	Mean Copies/Reaction
22Rv1 Pos	3	14.98	0.06	582422.25
h.blood Neg	3	ND	ND	ND
NTC	3	ND	ND	ND

FIGURE 9B

Experimental (in human blood gDNA)

1-Stage qPCR_pET28b-XMRV env Spiked Human Blood gDNA						
Sample (human blood gDNA)	pET28b XMRV env	N Rows	Mean Ct	Standard Deviation	Copies/ reaction	%Recovery of Copies/reaction
Neat	1E0	3	33.44	0.51	0.68	68.00
	1E1	3	30.12	0.07	7.62	76.20
	1E2	3	26.97	0.07	79.21	79.21
	1E3	3	23.82	0.13	825.75	82.58
	1E4	3	20.51	0.08	9641.13	96.41
	1E5	3	16.89	0.1	141843.44	141.84
	1E6	3	14.85	0.04	643411.76	64.34
	1E7	3	12.23	0.07	4494825.98	44.95

FIGURE 9C

XMRV pol2
Standards

1-Stage qPCR_pAZ3-emd pol2 Plasmid Standards			
pAZ3-emd pol2 Input Copies/Reaction	N Rows	Mean Ct	Standard Deviation
1E0	3	35.15	1 of 3 Detected
1E1	3	33.15	1.15
1E2	3	28.99	0.1
1E3	3	26.36	0.11
1E4	3	23.36	0.08
1E5	3	20.73	0.14
1E6	3	17.57	0.1

FIGURE 10A

Controls

1-Stage qPCR_XMRV pol2 Controls				
Controls	N Rows	Mean Ct	Standard Deviation	Mean Copies/Reaction
22Rv1 Pos	3	14.62	0.08	9341829.77
h.blood Neg	3	ND	ND	ND

FIGURE 10B

Experimental (in human blood gDNA)

1-Stage qPCR_pAZ3-emd pol2 Spiked Human Blood gDNA						
Sample (human blood gDNA)	pAZ3-emd pol2	N Rows	Mean Ct	Standard Deviation	Copies/ reaction	%Recovery of Copies/reaction
Neat	1E0	3	34.46	0.1	2.21	221.00
	1E1	3	34.13	0.91	3.3	33.00
	1E2	3	29.74	0.24	84.01	84.01
	1E3	3	26.69	0.06	866.55	86.66
	1E4	3	23.97	0.05	7011.98	70.12
	1E5	3	20.98	0.07	70064.46	70.06
	1E6	3	17.67	0.02	885681.45	88.57

FIGURE 10C

0-Stage vs. 1-Stage qPCR Protocols
XMRV-Gag
pUC57 XMRV gag_Standards

pUC57 XMRV gag_Standards				
pUC57 XMRV gag	Protocol Stage	N Rows	Mean Ct	Standard Deviation
1E0	0-stage	3	37.03	0.02
	1-stage	3	35.15	1.42
1E1	0-stage	3	33.76	0.16
	1-stage	3	30.83	0.1
1E2	0-stage	3	30.41	0.06
	1-stage	3	27.45	0.14
1E3	0-stage	3	26.32	0.1
	1-stage	3	23.32	0.08
1E4	0-stage	3	23.94	0.07
	1-stage	3	20.87	0.09
1E5	0-stage	3	20.13	0.02
	1-stage	3	17.13	0.06
1E6	0-stage	3	17.17	0.03
	1-stage	3	14.47	0.09
1E7	0-stage	3	15.66	0.05
	1-stage	3	12.71	0.02

FIGURE 11A

pUC57 XMRV gag spike-ins into 001 Human Blood gDNA

pUC57 XMRV gag spike-ins into 001 Human Blood gDNA					
pUC57 XMRV gag/001 gDNA	Protocol Stage	N Rows	Mean Ct	Standard Deviation	Detected Copies/reaction
1E0	0-stage	3	36.64	0.94	1.07
	1-stage	3	34.02	0.52	1.1
1E1	0-stage	3	32.96	0.18	13.97
	1-stage	3	29.53	0.69	27.73
1E2	0-stage	3	30.23	0.08	101.91
	1-stage	3	27.26	0.13	127.75
1E3	0-stage	3	26.5	0.08	1553.29
	1-stage	3	23.52	0.08	1826.01
001 Blood gDNA_Naive	0-stage	3	ND	ND	ND
	1-stage	3	ND	ND	ND

FIGURE 11B

EP 2 913 402 B1

XMRV-Env
pET28b XMRV env_Standards

pET28b XMRV env_Standards				
pET28b XMRV env	Protocol Stage	N Rows	Mean Ct	Standard Deviation
1E0	0-stage	3	36.08	0.17
	1-stage	3	32.9	0.63
1E1	0-stage	3	32.34	0.21
	1-stage	3	29.8	0.36
1E2	0-stage	3	28.82	0.16
	1-stage	3	26.59	0.11
1E3	0-stage	3	25.7	0.05
	1-stage	3	23.53	0.07
1E4	0-stage	3	22.4	0.05
	1-stage	3	20.22	0.04
1E5	0-stage	3	18.66	0.09
	1-stage	3	16.5	0.07
1E6	0-stage	3	16.56	0.15
	1-stage	3	14.14	0.07
1E7	0-stage	3	14.35	0.1
	1-stage	3	12.11	0.05

FIGURE 12A

pET28b XMRV env spike-ins into 001 Human Blood gDNA

pET28b XMRV env spike-ins into 001 Human Blood gDNA					
pET28b XMRV env /CA gDNA	Protocol Stage	N Rows	Mean Ct	Standard Deviation	Detected Copies/reaction
1E0	0-stage	3	36.42	0.74	0.51
	1-stage	3	34.71	0.47	0.23
1E1	0-stage	3	33.49	0.79	4.53
	1-stage	3	30.43	0.27	5.59
1E2	0-stage	3	29.88	0.1	56.34
	1-stage	3	27.08	0.04	68.29
1E3	0-stage	3	25.86	0.3	1078.16
	1-stage	3	23.16	0.08	1301.04
001 Blood gDNA_Naive	0-stage	3	ND	ND	ND
	1-stage	3	ND	ND	ND

FIGURE 12B

EP 2 913 402 B1

XMRV-Pol2

pAZ3-emd pol2_Standards

pAZ3-emd pol2_Standards				
pAZ3-emd pol2	Protocol Stage	N Rows	Mean Ct	Standard Deviation
1E0	0-stage	3	ND	ND
	1-stage	3	36.06	1 of 3 Detected
1E1	0-stage	3	36.01	0.3
	1-stage	3	33.09	0.15
1E2	0-stage	3	32.75	0.16
	1-stage	3	29.99	0.2
1E3	0-stage	3	29.79	0.13
	1-stage	3	26.56	0.04
1E4	0-stage	3	26.7	0.09
	1-stage	3	23.47	0.05
1E5	0-stage	3	23.93	0.06
	1-stage	3	20.57	0.09
1E6	0-stage	3	20.99	0.14
	1-stage	3	17.75	0.16

FIGURE 13A

pAZ3-emd pol2 spike-ins into 001 Human Blood gDNA

pAZ3-emd pol2 spike-ins into 001 Human Blood gDNA

pAZ3-emd pol2 /001 gDNA	Protocol Stage	N Rows	Mean Ct	Standard Deviation	Detected Copies/reaction
1E0	0-stage	3	ND	ND	ND
	1-stage	3	35.05	1 of 3 Detected	2.09
1E1	0-stage	3	36.74	0.73	5.44
	1-stage	3	33.63	1.22	7.38
1E2	0-stage	3	32.65	0.01	114.29
	1-stage	3	29.52	0.18	130.53
1E3	0-stage	3	30.06	0.21	847.39
	1-stage	3	27.18	0.08	746.81
001 Blood gDNA_Naive	0-stage	3	37.49	1 of 3 Detected	2.72
	1-stage	3	ND	ND	ND

FIGURE 13B

EP 2 913 402 B1

Sample	Spike-In Copy #	MLV					Env2				
		Mean Ct	Std Dev	N	% CV	Copy #	Mean Ct	Std Dev	N	% CV	Copy #
100 ng Frozen Fresh +	7980	26.05	0.065	3	0.25%	2660	25.05	0.119	3	0.47%	3568
100 ng FFPE +	4856	27.12	0.434	3	1.60%	1301	24.28	0.408	3	1.68%	6229
100 ng FFPE -	-----	27.48	1.060	3	3.86%	1141	25.39	0.676	3	2.66%	3029

FIGURE 14

Input Copies/reaction	Standards	N Rows	Mean Ct	Standard Deviation	%CV
1.00E+00	XMRV env	3	.	.	.
	XMRV gag	2	.	.	.
	XMRV pol2	3	.	.	.
1.00E+01	XMRV env	3	32.92	0.77	2.35
	XMRV gag	2	30.83	0.28	0.89
	XMRV pol2	3	33.63	0.69	2.06
1.00E+02	XMRV env	3	29.43	0.19	0.66
	XMRV gag	2	30.83	0.28	0.89
	XMRV pol2	3	30.04	0.05	0.17
1.00E+03	XMRV env	3	26.58	0.11	0.4
	XMRV gag	2	27.6	0.33	1.18
	XMRV pol2	3	26.61	0.2	0.75
1.00E+04	XMRV env	3	23.4	0.06	0.24
	XMRV gag	2	24.3	0.09	0.38
	XMRV pol2	3	23.41	0.07	0.31
1.00E+05	XMRV env	3	20.56	0.04	0.2
	XMRV gag	2	21.45	0.13	0.59
	XMRV pol2	3	20.86	0.24	1.13
22Rv1 gDNA_5 ng_Pos	XMRV env	3	19.28	0.09	0.47
	XMRV gag	3	20.39	0.1	0.47
	XMRV pol2	3	19.71	0.18	0.91
qPCR NTC	XMRV env	3	.	.	.
	XMRV gag	3	.	.	.
	XMRV pol2	3	.	.	.
RT NTC	XMRV env	3	.	.	.
	XMRV gag	2	.	.	.
	XMRV pol2	3	.	.	.

FIGURE 15




FIGURE 16A




FIGURE 16B

REFERENCES CITED IN THE DESCRIPTION

This list of references cited by the applicant is for the reader's convenience only. It does not form part of the European patent document. Even though great care has been taken in compiling the references, errors or omissions cannot be excluded and the EPO disclaims all liability in this regard.

Patent documents cited in the description

- WO 2010036986 A [0006] [0011] [0097]
- WO 2010075414 A [0007]
- WO 2010045002 A [0011]
- US 2009058512 W [0025]
- WO 9822489 A [0029]
- WO 9839352 A [0029]
- WO 9914226 A [0029]
- US 5143854 A [0045] [0046]
- WO 9015070 A [0045] [0046]
- WO 9210092 A [0045] [0046]
- US 5412087 A [0046]
- WO 9511995 A [0046]
- US 4965188 A [0056]
- WO 9001069 A [0056]
- US 5130238 A [0056]
- EP 329822 A [0056]
- US 5169766 A [0056]
- WO 8906700 A [0056]
- US 4458066 A [0065]
- US 0958512 W [0097]
- WO 2010148203 A [0113]
- US 61365297 B [0118]
- WO 61386941 A [0118]
- WO 61391360 A [0118]

Non-patent literature cited in the description

- **URISMAN et al.** *PLOS pathogens*, 2006, vol. 2 (3), e25 [0002]
- **DONG et al.** *PNAS*, 2007, vol. 104, 1655 [0002]
- **R. SCHLABERG et al.** *PNAS*, 2009 [0003]
- **V.C.LOMBARDI et al.** *Science*, 2009 [0003]
- **FISCHER et al.** *Journal of Clinical Virology*, 2008, vol. 43, 277-283 [0003]
- **HOHN et al.** *Retrovirology*, 2009, vol. 6, 92 [0003]
- **VAN KUPPEVELD et al.** *BMJ*, 2010, vol. 340, c1018 [0003]
- **FISCHER et al.** *Emerg Infect Dis.*, 2010 [0003]
- **DONG et al.** *PNAS*, 2007 [0004]
- **ENSERINK et al.** *Science*, 2010, vol. 329, 18-19 [0004]
- **PAPROTKA T.** *Science*, 2011, vol. 333, 97-101 [0005]
- **ROBINSON MJ. et al.** *Retrovirology*, 2010, vol. 7, 108 [0005]
- **TAI et al.** *Mol. Ther.*, 2005, vol. 12, 842-851, http://oba.od.nih.gov/oba/RAC/meetings/Jun2009/976_Aghi.pdf [0006]
- **GURUNATHAN S ; HABIB RE et al.** *Vaccine*, 2009, vol. 27, 1997-2015 [0006]
- **LOW A. ; OKEOMA CM et al.** *Virology*, 2009, vol. 385, 455-463 [0006]
- **URISMAN et al.** Identification of a novel Gammaretrovirus in prostate tumours of patients homozygous for R462Q RNASEL variant. *PLOS Pathogens*, 2006, vol. 2 (3), E25-1 [0008]
- **LOGG et al.** *J. Virol.*, 2001, vol. 75, 6989-6998 [0011]
- **SAKUMA et al.** *Virology*, 2009 [0016]
- **POWELL et al.** *J. Virol.*, 1999, vol. 73, 8813-8816 [0016]
- **G.B.BECK-ERIGESER.** *PNAS*, 2009 [0016]
- **J.D. FOX.** *J.Clin. Virol.*, 2007, vol. 40 (1), S15-S23 [0018]
- **ANASTASSOPOULOU et al.** *Journal of Virological Methods*, 2001, vol. 91, 67-74 [0018]
- **GARBESI.** *Nucl. Acids Res.*, 1993, vol. 21, 4159-65 [0030]
- **FUJIMORI.** *J. Amer. Chem. Soc.*, 1990, vol. 112, 7435 [0030]
- **URATA.** *Nucleic Acids Symposium Ser.*, 1993, vol. 29, 69-70 [0030]
- **KORNBERG ; BAKER.** *DNA Replication.* Freeman, 1992 [0030]
- **SAMBROOK et al.** *Molecular Cloning, A Laboratory Manual.* Cold Spring Harbor Press, 1989 [0052]
- *Nucleic Acid Hybridization.* **HAYMES et al.** *A Practical Approach.* IRL Press, 1985 [0052]
- **MANIATIS et al.** *Molecular Cloning: A Laboratory Manual.* 1989 [0053]
- *Short Protocols in Molecular Biology* [0053]
- **BARANY et al.** *Proc. Natl. Acad. Sci. USA*, 1991, vol. 88, 189-93 [0056]
- **LANDEGREN et al.** *Science*, 1988, vol. 241, 1077-80 [0056]
- **WALKER et al.** *Proc. Natl. Acad. Sci. USA*, 1992, vol. 89, 392-6 [0056]
- **NARANG et al.** *Methods in Enzymology*, 1979, vol. 68, 190 [0065]
- **BROWN et al.** *Methods in Enzymology*, 1979, vol. 68, 109 [0065]
- **BEAUCAGE et al.** *Tetrahedron Letters*, 1981, vol. 22, 1859 [0065]

- **E.C. KNOUF et al.** *J. Virol*, 2009, vol. 83, 78353-7356
[0091]