



(72) ANDERSSON, Birger, SE

(71) ANDERSSON, Birger, SE

(51) Int.Cl.<sup>6</sup> C12Q 1/18, G01N 33/68, G01N 33/58

(30) 1995/11/01 (9502409-7) SE

(54) **PROCEDE D'ANALYSE IN VITRO DE SUBSTANCES  
TOXIQUES ET ALLERGENES**

(54) **A PROCESS FOR IN VITRO ANALYSIS OF TOXIC AND  
ALLERGENIC SUBSTANCES**

(57) Procédé d'évaluation in vitro d'une substance potentiellement allergène ou pouvant irriter les tissus, selon lequel des cellules sanguines sont cultivées en présence de dilutions successives de ladite substance, la plus forte concentration de substance étant non toxique pour les cellules est diluée de manière successive, et la prolifération cellulaire et la présence de cytokines sont mesurées. Selon le présent procédé, la présence d'une ou plusieurs des cytokines d'alarme de classe 0 uniquement est une indication de lésions de tissus et d'effets chimiques toxiques, la présence d'une ou plusieurs cytokines d'alarme de classe 0 et éventuellement d'une ou plusieurs cytokines de classe IV, à l'exception de la néoptérine, est une indication d'une hypersensibilité retardée, telle que l'immunité cellulaire, l'allergie retardée et l'eczéma de contact et la présence d'une ou plusieurs cytokines d'alarme de classe 0, d'au moins la néoptérine et éventuellement d'une ou plusieurs cytokines de classe I est une indication d'hypersensibilité immédiate telle que l'asthme, le rhume des foies, l'urticaire et la rhinite. La présente invention concerne l'utilisation de la néoptérine et de l'interleukine-8 en tant que substances d'analyse pour distinguer entre les réactions immunes des types I et IV, ainsi qu'un kit d'analyse.

(57) The invention concerns a process for in vitro evaluation of a potentially allergenic or tissue irritating substances, characterized in that blood cells are cultivated in the presence of serial dilutions of the substance whereby the highest concentration of the substance being non toxic to the cells is serially diluted, cell proliferation is measured and the presence of cytokines is measured whereby: the presence of one or more of the alarm cytokines of class 0 only is an indication of tissue damage and chemical toxic effects; the presence of one or more alarm cytokines of class 0 and possibly one or more cytokines of class IV type but not neopterin is an indication of delayed type hypersensitivity such as cellular immunity, delayed allergy and contact eczema; and the presence of one or more alarm cytokines of class 0 and at least neopterin and possibly one or more cytokines of class I is an indication of immediate type hypersensitivity such as asthma, hay fever, urticaria and rhinitis. The invention also concerns use of neopterin and IL-8 as analytic substances for distinguishing between the immune reactions type I and IV and a kit for analysis.



**PCT**WORLD INTELLECTUAL PROPERTY ORGANIZATION  
International Bureau

## INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

<b>(51) International Patent Classification <sup>6</sup> :</b> <b>G01N 33/68, 33/50, C12Q 1/02, C12N 5/08</b>	<b>A1</b>	<b>(11) International Publication Number:</b> <b>WO 97/16732</b> <b>(43) International Publication Date:</b> 9 May 1997 (09.05.97)
<b>(21) International Application Number:</b> PCT/SE96/00902 <b>(22) International Filing Date:</b> 3 July 1996 (03.07.96) <b>(30) Priority Data:</b> 9502409-7 1 November 1995 (01.11.95) SE <b>(71)(72) Applicant and Inventor:</b> ANDERSSON, Birger [SE/SE]; Jungfrugatan 45, S-114 44 Stockholm (SE). <b>(74) Agents:</b> BERG, S., A. et al.; H. Albihns Patentbyrå AB, P.O. Box 3137, S-103 62 Stockholm (SE).		<b>(81) Designated States:</b> AL, AM, AT, AU, AZ, BB, BG, BR, BY, CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GE, HU, IL, IS, JP, KE, KG, KP, KR, KZ, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, TJ, TM, TR, TT, UA, UG, US, UZ, VN, ARIPO patent (KE, LS, MW, SD, SZ, UG), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG).  <b>Published</b> <i>With international search report.</i>
<b>(54) Title:</b> A PROCESS FOR IN VITRO ANALYSIS OF TOXIC AND ALLERGENIC SUBSTANCES <b>(57) Abstract</b> <p>The invention concerns a process for in vitro evaluation of a potentially allergenic or tissue irritating substances, characterized in that blood cells are cultivated in the presence of serial dilutions of the substance whereby the highest concentration of the substance being non toxic to the cells is serially diluted, cell proliferation is measured and the presence of cytokines is measured whereby: the presence of one or more of the alarm cytokines of class 0 only is an indication of tissue damage and chemical toxic effects; the presence of one or more alarm cytokines of class 0 and possibly one or more cytokines of class IV type but not neopterin is an indication of delayed type hypersensitivity such as cellular immunity, delayed allergy and contact eczema; and the presence of one or more alarm cytokines of class 0 and at least neopterin and possibly one or more cytokines of class I is an indication of immediate type hypersensitivity such as asthma, hay fever, urticaria and rhinitis. The invention also concerns use of neopterin and IL-8 as analytic substances for distinguishing between the immune reactions type I and IV and a kit for analysis.</p>		

A process for in vitro analysis of toxic and allergenic substances

5 The present invention concerns a process for in vitro  
analysis of toxic and allergenic substances. Substances  
intended for use as pharmaceuticals, food additives, cosmetic  
or hygienic products, industrial chemicals and other sub-  
stances are analysed for adverse reactions. Such reactions  
may be allergic reactions, skin irritating effects and toxic  
10 effects.

Today such tests are performed in vivo on animals. Animals  
are expensive to raise and keep. Further, some tests may  
involve suffering of the animals.

15 The present invention provides a process for in vitro  
analysis of the adverse effects of substances. The analysis  
of the substances is performed on blood cells from warm  
blooded animals. No animals are involved since the substances  
are tested on human cells i.e. cells from the same species as  
20 they are intended for. Some substances may not have the same  
effect on different species, and tests performed on animals  
may give false results. Therefore the test is preferably  
performed with human blood.

25 Thus, the method will be offered as an alternative to animal  
tests for food additives, cosmetic or hygienic products,  
industrial chemical, drugs and other substances where adverse  
reactions are to be avoided.

30 In vitro evaluation of toxic effects on the immune system has  
been tested on splenocytes from mice ("In vitro evaluation of  
drug-induced toxic effects on the immune system as assessed  
by proliferative assays and cytokine production", M. Pallardy  
35 et. al. Eur. Cytokine Net., Vol. 2 No 3, May-June 1991, pp.  
201-206). The method was poorly effective detecting molecules  
inducing autoimmunity and hypersensitivity.

It has now turned out that according to the present invention it is possible to analyse adverse tissue reactions to foreign substances in vitro by cultivating them with human whole blood, prevented from coagulating and analysing the cytokine produced by the human cells.

Before a substance has been in use and before the target group has been exposed the present method offers the possibility to predict whether a risk for adverse reactions of inflammatory or allergic type will occur.

The invention provides an interpretation key for evaluating the adverse tissue reactions to foreign substances.

Summary of the invention: The method is defined in claim 1 and is based on cultivation of blood cells in vitro. The substance to be evaluated is added to the blood culture. If the reactive cells of the blood recognises the test material as foreign or damaging a production of inflammatory signal substances, cytokines occurs. The blood contains different white blood cells with different functions in the response to foreign substances. The different types of inflammatory cells secrete different repertoires of cytokines. The patterns, (profiles) of the secreted cytokines indicate what compartment of the immune system that has been activated. Thereby is it also possible to predict what kind of adverse reaction to be expected from a substance if administered to humans or experimental animals. Some substances will induce a typical cytokine profile driving the immune system towards IgE mediated allergy, thus indicating a potential of those substances to evoke asthma, hay fever and urticaria. Other substances will induce a cytokine pattern driving the immune system to cellular immunity thus having the potential to result in contact excema as an adverse reaction. Finally, a chemical or toxic effect will result in a cytokine pattern typical of alarm reactions seen in tissue damage and operations. Thus, the concept that the present method is based on is to use an interpretation of the cytokine pattern that a

substance can induce in blood in vitro cells to predict whether the substance has the potential to cause allergy of different types or is generally hazardous to induce tissue damaging or toxic effects.

We classify adverse tissue reactions to foreign substances as follows.

Class 0: Tissue damage. Unspecific chemical, mechanical or physical damage.

Class I: Allergic immune reaction type I. This is also named immediate type hypersensitivity and is mediated by IgE antibodies produced specifically to the foreign substance. An acute inflammatory reaction is produced, histamine is often produced and examples of symptoms are asthma hay fever, urticaria and rhinitis. This type requires that the substance can induce a fully mature immune response where all components of the immune system participates. This is also considered as a final step immune reaction.

Class IV: Inflammatory immune reaction type IV. This is also named delayed type hypersensitivity and is mediated by sensitised T lymphocytes type TH-1, considered as a first step type immune reaction. Allergic contact dermatitis is an example of this type of reaction.

Class 0 is our own classification of unspecific damage to the immune system. For Class I and Class IV see "Immunology" by Ivan Roitt, Jonathan Brostoff and David Male, Gower Medical Publishing London - New York, 1989, page 19.1-19.20 and 22.1-22.10, which is incorporated as a reference.

Three main types of cytokine profiles have been identified in the present work.

Class 0 type: Secretion of alarm cytokines indicating damage to connective tissue, fibroblasts, endothelial cells,

epithelial cells and unspecific inflammatory white blood cells. Members of this group are IL-1, IL-6, IL-12 and TNF.

Class I type: Secretion of cytokines of an immune response type from lymphocytes and inflammatory cells. These includes cytokines from Class 0 and in addition IFN-gamma, Neopterin, IL-2 and IL-10. Theoretically as known from animal in vivo studies also IL-4 and IL-5 should be produced here, but these substances are notoriously difficult to determine and are therefore not routinely included in our test protocols.

Class IV type: Secretion of the cytokines of Class 0 unspecific type and in addition IL-2, IL-8, IL-10 and IFN-gamma.

A finding that could not predicted was that the complete final type immune response Class I drives to Neopterin production whereas the Class IV primary type immune response does not stimulate the immune system as long as to Neopterin production.

In the three groups of cytokines listed above only preferably selected representative members are listed. The groups may contain other cytokines to be used and the patent is claimed for the use of analysis of all substances listed in The Cytokine Facts Book, Eds Callard R.E. and Gearing, A.J.H., Academic Press 1994 hereby incorporated as a reference and future updated issues to predict the grade of maturity of the response that a substance has the potential to evoke. At present 50 substances.

#### **Interleukins**

IL-1, IL-2, IL-3, IL-4, IL-5, IL-6, IL-7, IL-8, IL-9, IL-10, IL-11, IL-12, IL-13, IL-14, IL-15

#### **Other Cytokines (in alphabetical order)**

BDNF, CNTF, EGF, Epo, FGF, G-CSF, GM-CSF, I-309/TCA-3,  $\gamma$ IP-10, IFN $\alpha$ , IFN $\beta$ , IFN $\gamma$ , LIF, LT(TNF $\beta$ ), MCP-1,2 and 3,

M-CSF, MIF, MIP-1 $\alpha$ , MIP-1 $\beta$ , MIP-2, NGF, NT-3, NT-4, OSM, PBP, PBSF, PDGF, PF-4, RANTES, SCF, TGF $\alpha$ , TGF $\beta$ , TNF $\alpha$ , Tpo, VEGF

According to the invention the presence of inflammatory immune reaction type IV is preferably analysed by using IL-8. It has turned out that IL-8 is developed and is present in raised levels during type IV reactions.

Neopterin being present in higher amounts when type I reaction is at hand, IL-8 and/or neopterin can especially be used to distinguish between the inflammatory immune reaction type IV and the allergic immune reaction type I.

Based on the outcome of the in vitro test, the cytokine pattern induced by a the substance under investigation, a prediction will be made regarding the potential of the substance to cause adverse reaction if humans or animals are exposed to it. Whole blood, preferably venous blood is used. The blood is prevented from coagulating. This can be done by the use of endotoxin free heparin tubes. Other methods may also be used such as addition of EDTA or citrate. Alternatively the blood can be shaken with glass beads.

The leukocyte counts in undiluted blood ranges between  $3 \times 10^9$  and  $7 \times 10^9$  per litre. Around 30% of these are lymphocytes and 5-10% are monocytes the rest being polymorfonuclear cells. Total leukocyte numbers and differential count is determined in a Coulter STKS cell haematology analyser before further procedure.

Whole blood is diluted serially from 1:1-1:1000, preferably 1:2-1:100 in a culture medium, preferably a tissue culture medium such as RPMI 1640 supplemented with glutamine. Normally 1,5-5mM, preferably 2mM is used. The blood is most often diluted 1:5, 1:10 1:20 or continued to 1:40, 1:80, 1:100.

Antibiotics are preferably added to prevent the growing of unwanted microorganisms. Penicillin and streptomycin may be added in concentrations of 25-100U/ml, preferably 50U/ml.

5 Mitogens or substances with known effects on the immune system can be used as positive controls to give better results. Any mitogen as mentioned in Daniel P. Stites: Clinical Laboratory Methods for Detection of Antigens & Antibodies in. Basic and Clinical Immunology, Lange Medical Publications, Los Altos California, 1984 may be used. This  
10 reference is incorporated in this description by reference. Phytohemagglutinin (PHA-L) is preferably used and dissolved in the medium to a final concentration of in the wells of 250 µg per ml.

15 The test can be performed in any suitable vessel e.g. in tubes or preferably in the wells of microtiter plates.

20 For triggering cytokine production in vitro the cell concentration is of critical importance. If the test is performed only at one cell concentration, then only some of the cytokines of interest will be optimally detected. If a series of cell concentration is used the complete picture of the cytokine pattern predicted above in " Summary of the  
25 invention" is revealed.

30 For substances to be tested for effects in the system different concentrations are added to the blood. The highest concentration of the substance being non toxic to the cells indicated by inspection in microscope in presence of vital stain e.g. trypan blue is used as start concentration. Serial dilution of the substance is then performed down to a dilution of at least 1:1 000 and preferably down to 1:10 000, 1:100 000 or 1:1 000 000.

35 In cases where the substance is diluted in RPMI 1640 the appropriate control is medium only. When other diluents are



used, controls of the appropriate diluent at corresponding concentrations are used.

The vessels are incubated at 25-40°C, preferably at 37°C at 50-95%, preferably 90% humidity and preferably in the presence of with 2-15% CO<sub>2</sub>, preferably 5%. An appropriate number of tubes or microtiter plates are prepared and are taken out of the incubator for testing at the intervals. Standard time intervals are 1 hour, up to 10 days; preferably 1 hour 1 day, 2 days and 4-6 days.

Proliferation assay can be made as follows. One set of plates is used to estimate the proliferative activity of the donor cells spontaneously ( =medium ) or substance dissolvent controls ), in response to mitogens and in response to samples of test substances at various concentrations This can be done as follows or by any other suitable method. At time intervals plates are tested, and to each well in the plate is added labeled nucleic acid precursor such as <sup>3</sup>H-thymidine dissolved in the medium. The plates are incubated for additional 4 hours at 37°C for incorporation of labeled nucleic acid precursor. The cells are then collected using a Skatron cell harvester (Skatron, Lier, Norway) and Titertek paper filters or any other suitable method. The result of this procedure is that the DNA of the cells is captured in the filter paper, whereas the <sup>3</sup>H-thymidin not incorporated in the DNA passes through the filter and goes with the waste washing water. The radioactivity in the filter papers is thus proportional to the DNA synthesis that occurred during the 4 hours preceding the harvesting procedure. The radioactivity of the filter papers is determined in a Packard liquid scintillation counter, and thymidine incorporation is expressed as counts per minute ( cpm ) The radioactivity or any other label of the precursor can be measured in DNA harvested by any other mechanical or chemical way and measured by any suitable isotopic or chemical instrument.

Cytokine assays: Another set of plates is prepared with the same design as in "Culture protocol" above. Also here, plates are removed from the incubator at various time intervals, and cytokines released from the cultured cells are measured in the supernatants obtained after centrifugation of the plates at 650 g for 10 minutes. The supernatants can either be tested immediately or stored at -20°C until tested. For cytokine determinations diagnostic kits of various origins are used. Patent is claimed for this procedure using any available or newly constructed test using cytokine quantitation including bioassays, immunoassays or chemical assays or other.

The instructions of the manufacturers are followed. All these tests are enzyme immunoassays ( EIA:s ). The principle is that microtiter plate wells are labeled with a cytokine specific capture antibody. If cytokine is present in the samples added to the wells it will be captured to the bottom of the well. To quantitate how much cytokine that has been captured a second antibody, labeled with an enzyme is added to the wells. The reaction is thereafter measured as colour developed after addition of a substrate for the enzyme. The values for the test wells is compared with a standard curves obtained from a series of known amounts of cytokine. The values for control wells with no substance added to the cells are used as background levels for the tests. Positive indications are recorded when values above 2 SD of the variation are obtained.

The invention also concerns a kit comprising one or more reagents recognizing the alarm cytokines of class zero, and/or the cytokines of class I type, and/or the cytokines of class IV. Examples of such cytokines are mentioned in the last two paragraphs on page 4 and the first paragraph on page 5. Reagents reacting for the presence of neopterin and IL-8 are preferred.

These reagents may be antibodies or any other reagents that are sensitive for these substances. The reagents may be supported by a carrier such as a strip, titer plates, micro-titer plates, Eliza plates, test tubes etc. The carriers may be of different sizes. The carrier material may be any solid or semi-solid material that does not interfere with the reaction between the reagent and the cytokine. The carrier can take on a variety of shapes and compositions, including microparticles, beads, porous and impermeable strips and membranes, the interior surface of reaction vessels such as test tubes and microtiter plates, and the like. Microtiter plates and beads may be of plastic, such as styrene or acryl-polymers or glass. Nitrocellulose can be used, preferably in the form of filters, strips or discs. Means for attaching a desired reaction partner to a selected solid support will be a matter of routine skill to the worker in the field. It is also possible to use flow cytometer.

The invention will now be better described with the following non limiting examples.

#### Example

Blood donors and blood collection: Five healthy blood donors from both sexes ranging in age between 25 and 55 years are used. Blood is collected by venous puncture into endotoxin free heparin tubes ( Coatech, Ljungby, Sweden ) that contains 120 IU heparin and are filled with 4 ml of blood. None of the blood donors are allowed to have a history of allergy or atopic excema or any indication of ongoing disease of any kind.

Cell analysis: Total and differential white blood cell counts are measured with Coulter STKS equipment. The five donors had total white blood cell counts of 9,8; 7,1; 9,8; 6,9; 5,8 · 10<sup>9</sup>/l.

Culture protocoll: The whole blood is diluted serially 1:2, 1:3, 1:20, 1:4, and 1:10 in tissue culture medium, RPMI 1640

supplemented with 2 mM glutamine, 50U/ml of penicillin and 50Ug/ml of streptomycin. 100 ul of blood is added to the wells of Nunclon microtiter plates. The mitogen phytohemagglutinin ( PHA-L ) ( L-4144, Sigma St Louis USA ) is dissolved in RPMI 1640 and 50 ul of the solution is added to the wells so that a final amount of 25 µg per well is achieved.

For substances to be tested for effects in the system different concentration are added to the blood. The highest concentration of the substance being non toxic to the cells indicated by inspection in microscope in presence of vital stain e.g. trypan blue is used as start concentration. The following substances and concentrations were used:

0 = Triton x 100 in concentrations 1 : 2000, 1 : 4000,

1 : 8000, 1 : 16000, 1 : 32000, 1 : 64000, 1 : 128000

IV = NiCl<sub>2</sub> 100 µg, 10 µg, 1 µg, 0,1 µg, 0,01 and 0,001 µg per ml

I = mite extract (1 : 10, 1 : 20, 1 : 40, 1 : 80

Pos = PHA 25 µg per well (200 µl)

In cases where the substance is diluted in RPMI 1640 the appropriate control is medium only. When other diluents are used, controls of the appropriate diluent at corresponding concentrations are used.

Incubation: The microtiter plates are incubated at 37°C at 90% humidity and with 5% CO<sub>2</sub>. An appropriate number of identical plates are prepared and are taken out of the incubator for testing at intervals of, 1 day, 2 days and 6 days.

Proliferation assay: One set of plates is used to estimate the proliferative activity of the donor cells spontaneously (=medium or substance dissolvent controls), in response to PHA, and in response to samples of test substances at various concentrations. The highest concentration of substance non toxic to the cells indicated by inspection in microscope in presence of vital stain (e.g. trypan blue) is used as start

concentration. Serial dilutions of the sub-stance is then performed as stated above.

To each well in the plate was added 1 uCi of  $^3\text{H}$ -thymidine (Amersham International, Amersham UK; 50 Ci/ mM) in 50 ul of RPMI 1640. The plates were incubated for addional 4 hours at 37°C. The cells were then collected using a Skatron cell harvester (Skatron, Lier, Norway) and Titertek paper filters. The result of his procedure is that the DNA of the cells is stuck in the filter paper, wheras the  $^3\text{H}$ -thymidin not incorporated in the DNA passes through the filter and goes with the waste washing water. The radioactivity in the filter papers is thus proportional to the DNA synthesis that occured during the 4 hours preceding the harvesting procedure. The radiactivity of the filter papers is determined in a Packard liquid scintillation counter, and thymidine incorporation was expressed as counts per minute (cpm).

Cytokine assays: Another set of plates is prepared with the same design as in "Culture protocoll " above. Cytokines released from the cultured cells are measured in the supernatants obtained after centrifugation of the plates at 650 g for 10 minutes. The supernatants were either tested immediately or stored at -20°C until tested.

For information only is here mentioned that in the tests performed the following commercial kits and manufacures were used.

IL-1beta: Immunotech, Chromgenix

IL-2: Immunotech, Chromgenix

IL-4: R&D Systems

Neopterin: Henning Berlin

IL-6: Immunotech, Chromgenix

IL-8: Assay Res. Inc

Interferon-gamma: Genzyme

sCD8: T Cell Diagnostics

sIL-2R Immunotech

TNF: Medgenix

IL-10: Medgenix

5 The instructions of the manufacturers were followed. All  
these tests are enzyme immunoassays ( EIA:s ). The principle  
is that microtiter plate wells are labeled with a cytokine  
specific capture antibody. If cytokine is present in the  
samples added to the wells it will be captured to the bottom  
10 of the well. To quantitate how much cytokine that has been  
captured a second antibody, labeled with an enzyme is added  
to the wells. The reaction is thereafter measured as colour  
developed after addition of a substrate for the enzyme. The  
values for the test wells is compared with a standard curve  
15 obtained by adding a series of known amounts of cytokine. The  
values obtained for control wells with no substance added are  
used as background values for cytokine levels and prolifera-  
tion as measured by <sup>3</sup>H-thymidin incorporation. Positive  
indications are recorded when values above 2 SD of the  
variation of the controls are obtained.

20 The results from testing three substances with previously  
from clinical experience known properties, but not previously  
tested with the method are shown:

- 25 0. Triton X-100 This is a detergent with tissue damaging  
effects. Dissolves cell membranes.
- I. Dust mite antigen ( Soluprick, ALK Sverige AB). This  
preparation is made from an extract of dust mites and  
is used clinically to test for immediate type allergy  
30 by intracutaneous injection.
- IV. Nickel chloride ( Sigma St Louis, USA ). This metal is  
wellknown as a cause of delayed type allergy ( contact  
excema ) in humans and experimental animals.
- Plus. The positive indicator of responsiveness of the system  
35 the calibrator was the mitogen PHA ( Sigma, St louis,  
USA ).

In vito performance test - Table I

Values for 11 cytokines and <sup>3</sup>H-Thymidine Proliferation. Obtained after stimulation of blood from 5 normal humans. All values are related to medium control given the value 100.

	IL-1	TNF	IL-6	IL-8	IL-2	STL-2R	SCD8	IL-10	IL-4	IFN $\gamma$	Neop- terin	Prolif. Index
O	630	400	460	200	85	115	95	95	85	120	110	200
IV	1165	225	660	380	445	130	110	800	80	1220	90	1340
I	320	240	250	100	900	120	85	640	900	800	560	120
Plus	690	420	1620	780	750	1820	940	1250	820	800	440	2600

The highest value as compared to medium control obtained at any time or concentration for each substance is used as indicator value in the table and related to medium control put to value 100 index.

5

Results: In summary the results of a performance example are presented below and the stimulatory effects were quite different for the different compounds.

10 0. Triton X-100 induced proliferation. Cytokines of the unspecific tissue toxic group were induced i.e. IL-1, IL-6 and TNF. This is a positive tissue irritant.

I. Dust mite: Induced no proliferative response.

15 Substances secreted into supernatants were the unspecific inflammatory IL-1, TNF and IL-6 and of immune lymphocyte class the IL-2, IL-4, IL-10 and IFN-gamma. Furthermore the Neopterin typical for induction of specific immune responses of a late mature stage. There is a potential for antibody and  
20 IgE production.

IV. Nickel chloride: Induced a low proliferative response. The cytokines induced were the unspecific inflammatory mediator IL-1, IL-6, TNF and IL-8. Of the immune lymphocyte  
25 class IFN-gamma, IL-2 and IL-10 were induced, but notably no Neopterin.

Plus. PHA: Induced a strong proliferative response indicated by <sup>3</sup>H-thymidine incorporation and also cell types of cytokines were induced. The substances secreted in the supernatants included: Cytokines of primary inflammatory type as  
30 IL-1, IL-6, TNF and IL-8. Cytokines of the types typical for specific immune responses both early and late as neopterin, IFN-gamma IL-2 and IL-4 and IL-10 typical of both early and late responses were induced.



## Claims

1. A process for in vitro evaluation of a potentially allergenic or tissue irritating substances, characterized in that blood cells are cultivated in the presence of serial dilutions of the substance whereby the highest concentration of the substance being non toxic to the cells is serially diluted, cell proliferation is measured and the presence of cytokines is measured whereby :

the presence of one or more of the alarm cytokines of class 0 only is an indication of tissue damage and chemical toxic effects,

the presence of one or more alarm cytokines of class 0 and possibly one or more cytokines of class IV type but not neopterin is an indication of delayed type hypersensitivity such as cellular immunity, delayed allergy and contact eczema and

the presence of one or more alarm cytokines of class 0 and at least neopterin and possibly one or more cytokines of class I is an indication of immediate type hypersensitivity such as asthma, hay fever, urticaria and rhinitis.

2. A process according to claims 1, characterized in that the cytokines of the alarm group are selected from cytokines indicating damage to connective tissue, fibroblasts, endothelial cells, epithelial cells and unspecific inflammatory white blood cells such as IL 1, IL-12, IL 6, TNF;

the cytokines of class IV type are selected from cytokines connected with cell mediated T-cells immunity such as IL-2, IL8, IL-10. IFN-gamma, IL-4, IL-5, and soluble products such as sCD8 and sIL-2R, but not neopterin;

the cytokines of class I type are selected among cytokines of an immune response type from T and B-lymphocytes and in-

flammatory cells such as IFN-gamma, IL-2, IL-10 , IL-4, IL-5 and soluble products such as sCD8 and sIL-2R and neopterin.

3. A process according to claim 1 and/or 2, characterized in  
5 that the substance to be tested is presented in the highest concentration being non toxic and in a serial dilution thereof at least 1:1000.

4. A process according to claim 3, characterized in that the  
10 serial dilution is performed down to a dilution of at least 1:10 000 to 1: 1 000 000.

5. A process according to claim 1 and/or 2, characterized in  
15 that the blood is diluted serially from 1:1 to 1:1000, preferably from 1:2 to 1:100.

6. A process according to claim 1 characterized in that  
whole blood is treated to not coagulate, diluted 1:2-1:100 in  
tissue culture medium preferably supplemented with glutamine  
20 and/or an antibioticum, a marked nucleic acid precursor is added,

the substance to be tested is serially diluted down to a  
dilution of 1: 1000 to 1:1 000 000 from the lowest con-  
25 centration of the substance being non toxic to the cells preferably as indicted by inspection in microscope in presence of vital stain and added in different concentrations and incubated at 25C-40C, and tested at time intervals from 1 hour to 6 days by

30 certifying that proliferation has taken place preferably by collecting the cells on filter paper or harvesting nucleic acid by any other mechanical or chemical method and analysing the labeled nucleic acid precursor and

35 analysing cytokines by estimating the presence of

one or more of the alarm cytokines of class 0 only is an indication of tissue damage and chemical toxic effects,

one or more alarm cytokines of class 0 and one or more  
5 cytokines of class IV type but not neopterin is an indication of delayed type hypersensitivity such as cellular immunity, delayed allergy and contact eczema and

one or more alarm cytokines of class 0 and at least neopterin  
10 and possibly one or more cytokines of class I is an indication of immediate type hypersensitivity such as asthma, hay fever and urticaria.

7. A process for in vitro evaluation of an potentially  
15 allergenic or tissue irritating substance characterized in that blood cells are cultivated in the presence of serial dilutions of the substance and in the presence of neopterin is analysed in a manner known per se whereby the presence of neopterin is an indication of immediate type hypersensitivity  
20 such as asthma, hay fever and urticaria.

8. Use of IL-8 and/or neopterin for analysing the presence of the inflammatory immune reaction type IV and/or the allergic immune reaction type I.  
25

9. A reagent kit for use in the process according to anyone of claims 1-7, characterized in that it comprises one or more reagents recognizing cytokines.

10. A reagent kit according to claim 9, characterized in that it comprises a reagent such as an antibody, recognizing IL-8 and/or neopterin.  
30