

(51) International Patent Classification:
A61K 35/16 (2006.01)(21) International Application Number:
PCT/US20 14/037496(22) International Filing Date:
9 May 2014 (09.05.2014)

(25) Filing Language: English

(26) Publication Language: English

(30) Priority Data:
61/821,3 19 9 May 2013 (09.05.2013) US
61/893,444 21 October 2013 (21.10.2013) US(71) Applicant: **ADVANCED NEUROREGENERATIVE THERAPIES, LLC** [US/US]; 10301 Hagen Ranch Road, Suite 600, Boynton Beach, FL 33437 (US).(72) Inventor: **MAHARAJ, Dipnarine**; 10301 Hagen Ranch Road, Suite 600, Boynton Beach, FL 33437 (US).(74) Agent: **LANDER, Ferris, H.**; Ferris H. Lander Inc., P.O. Box 31166, Palm Beach Gardens, FL 33420 (US).

(81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IR, IS, JP, KE, KG, KN, KP, KR, KZ, LA, LC, LK, LR, LS, LT, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SA, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.

(84) Designated States (unless otherwise indicated, for every kind of regional protection available): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, KM, ML, MR, NE, SN, TD, TG).

[Continued on nextpage]

(54) Title: METHOD OF REPAIRING AGE AND DISEASE IMMUNE DYSFUNCTION AND CELLULAR SENESENCE WITH LYMPHOID STEM CELLS AND THEN RE-APPLYING THOSE FOR THERAPEUTIC USE




FIGURE 2

(57) Abstract: This invention relates to methods for treatment of diseases of ageing including immunosenescence, immune dysfunction, inflammation and impairment of early lymphoid lineage differentiation. The invention more specifically relates to the use of granulocyte colony stimulating factors to assist in stem cell mobilization, optionally in combination with the application of a method of delivering precise magnetic field patterns which agree with the body's own natural magnetic field patterns, and further in combination with re-infusion of previously collected autologous cells and/or plasma, optionally including allogeneic (healthy donor) cells and blood plasma.



Published:

- *without international search report and to be republished
upon receipt of that report (Rule 48.2(g))*

**METHOD OF REPAIRING AGE AND DISEASE IMMUNE DYSFUNCTION AND
CELLULAR SENESENCE WITH LYMPHOID STEM CELLS AND THEN RE-APPLYING
THOSE FOR THERAPEUTIC USE**

CROSS-REFERENCE TO RELATED APPLICATIONS

[0001] This application claims benefit of the priority of U.S. Provisional Application Serial No. 61/821,319, filed on May 9, 2013, and of U.S. Provisional Application Serial No. 61/893,444, filed on October 21, 2013, the contents of each of which is incorporated by reference herein in its entirety.

FIELD OF THE INVENTION

[0002] This invention relates to methods for treatment of diseases of ageing including immunosenescence, immune dysfunction, inflammation and impairment of early lymphoid lineage differentiation. The invention more specifically relates to the use of granulocyte colony stimulating factors to assist in stem cell mobilization, optionally in combination with the application of a method of delivering precise magnetic field patterns which agree with the body's own natural magnetic field patterns, and further in combination with re-infusion of stem cell containing compositions, such as previously collected autologous cells and/or plasma, optionally including allogeneic (healthy donor) cells and blood plasma.

BACKGROUND OF THE INVENTION

[0003] The phenomenon which manifests as growth arrest after a period of apparently normal cell proliferation is known as Replicative Senescence (RS). Replicative Senescence is seen in a) cells from adults of all ages b) embryonic tissues, and c) animals.

[0004] Aging is associated with alterations of the immune system including impairments in innate immunity, T-lymphopoiesis and B-lymphopoiesis and these

impairments contribute to immunosenescence and immune dysfunction in affected individuals. Multipotent hematopoietic stem cells (HSCs) aging contributes to impairments in early lymphoid lineage differentiation. Immunosenescence with immune dysfunction and increased inflammation is a primary cause of aging and diseases such as anemia, chronic diseases, autoimmune disorders, cancer, cardiovascular diseases, infection, metabolic diseases, neurodegenerative diseases, protein energy malnutrition and frailty.

[0005] Oscillating magnetic fields have been used for years in the course of administering physical therapy to clinic patients suffering from bone fractures. These devices are typically called bone growth stimulators. Bone growth occurs as a result of stem cell stimulation, activation and differentiation. These device signals, however, are a series of pulses or oscillating waves, which have symmetry typical of electronic-generated signals (see Figure 1 "Common electronic-generated signals"). More recently, researchers have discovered that the body emits its own complex electromagnetic field pattern. Unique patterns are associated with immunosenescence and immune dysfunction, stress, or disease. By capturing these abnormal patterns, re-storing and re-admitting these patterns to the target patient, researchers theorize that the normal "healing process" may be restored more effectively, as the patterns would be natural biologic patterns.

[0006] What is unique about the instantly disclosed method described is the confluence of these unique processes to promulgate a therapeutic modality.

[0007] The number of Cumulative population doublings (CPDs) cells undergo in culture varies considerably between cell types and species. Early results suggested a relation between CPDs cells could endure and the longevity of the species from which the cells were derived, e.g. cells from the Galapagos tortoise, which can live over a century, divide about 110 times while mouse cells divide roughly 15 times. Cells taken from patients with progeroid syndromes such as Werner syndrome (WS)—exhibit far fewer CPDs than normal

humans. Certain "immortal" cell lines can divide indefinitely without reaching RS, e.g. embryonic germ cells and most cell lines derived from tumors, such as HeLa cells.

[0008] Biomarkers of cell senescence include:

- 1) Growth arrest - Senescent cells are growth arrested in the transition from phase G1 to phase S of the cell cycle. The growth arrest in RS is irreversible in the sense that growth factors cannot stimulate the cells to divide even though senescent cells can remain metabolically active for long periods of time;
- 2) Cellular morphology - Senescent cells are bigger and a senescent population has more diverse morphotypes than cells at earlier CPDs (Note Figure 12 which shows Normal human fibroblasts (left) and fibroblasts showing a senescent morphology (three cells on the right). Notice the common elongated morphology of senescent cells.
- 3) Senescence-associated β -galactosidase (SA β -gal) activity - In vitro and in vivo, the percentage of cells positive for SA β -gal increases with, respectively, CPDs and age. In immortal cell lines, such as HeLa tumor cells, the percentage of cells positive for SA β -gal does not correlate with CPDs. The increase in SA β -gal also correlates with the appearance of the senescent morphotypes;
- 4) Polyploid Increase - the percentage of polyploid cells—i.e., cells with three or more copies of chromosomes—has been shown to increase. Deletions in the mitochondrial DNA (mtDNA) have also been observed both during RS and during aging in vivo, at least in some cells;
- 5) Change in Gene Expression Levels - The expression levels of several genes change during in vitro cellular aging. One important type of gene

overexpressed in senescent cells are inflammatory regulators like interleukin 6 (IL6); proinflammatory proteins secreted by senescent cells in driving senescence, which may lead to positive feedback loops and to senescence induction in normal cells near senescent cells;

- 6) Metalloproteinase and Heat Shock Protein Production - Senescent cells also display an increased activity of metalloproteinases which degrade the extracellular matrix and a decreased ability to express heat shock proteins;
- 7) Telomere shortening - the primary cause of RS in human fibroblasts which have a major role in aging.

SUMMARY OF THE INVENTION

[0009] The phenomenon which manifests as growth arrest after a period of apparently normal cell proliferation is known as Replicative Senescence (RS). Replicative Senescence is seen in a) cells from adults of all ages b) embryonic tissues, and c) animals. The instant invention discloses treatment modalities for treating a number of maladies which result from the aging process.

[0010] Aging is associated with alterations of the immune system including impairments in innate immunity, T-lymphopoiesis and B-lymphopoiesis and these

impairments contribute to immunosenescence and immune dysfunction in affected individuals. Multipotent hematopoietic stem cells (HSCs) aging contributes to impairments in early lymphoid lineage differentiation. Immunosenescence with immune dysfunction and increased inflammation is a primary cause of aging and diseases such as anemia, chronic diseases, autoimmune disorders, cancer, cardiovascular diseases, infection, metabolic diseases, neurodegenerative diseases, protein energy malnutrition and frailty.

[0011] Oscillating magnetic fields have been used for years in the course of administering physical therapy to clinic patients suffering from bone fractures. These devices are typically called bone growth stimulators. Bone growth occurs as a result of stem cell stimulation, activation and differentiation. These device signals, however, are a series of pulses or oscillating waves, which have symmetry typical of electronic-generated signals (see Figure 1 "Common electronic-generated signals"). More recently, researchers have discovered that the body emits its own complex electromagnetic field pattern. Unique patterns are associated with immunosenescence and immune dysfunction, stress, or disease. By capturing these abnormal patterns, re-storing and re-admitting these patterns to the target patient, researchers theorize that the normal "healing process" may be restored more effectively, as the patterns would be natural biologic patterns. U.S. Patent 7,361,136 to Parker describes a method of treatment utilizing such a device, and is incorporated by reference herein in its entirety.

[0012] The instantly disclosed method describes a therapeutic modality that represents a confluence of one or more of these treatments.

[0013] Accordingly, it is a primary objective of the instant invention to treat diseases of ageing, diseases of ageing with immunosenescence and immune dysfunction and inflammation, and impairments in early lymphoid lineage differentiation by use of a stem cell mobilization agents, G-CSF, Granulocyte Colony Stimulating Factor) in combination with one or more of:

[0014] collecting autologous stem cells and plasma using a cell collection device;

[0015] re-infusing the previously collected autologous cells and or blood plasma;

[0016] re-infusing the previously collected autologous cells together with allogeneic (healthy donor) cells and or blood plasma; and

[0017] re-infusing allogeneic (healthy donor) cells and or blood plasma .

[0018] Treating diseases of ageing and diseases of ageing with immunosenescence and immune dysfunction and inflammation by stem cell activation using a method of delivering precise magnetic field patterns which agree with the body's own natural magnetic field patterns, and do so with an instrument capable of routine clinical therapy use, in combination with stem cell mobilization agents, G-CSF, Granulocyte Colony Stimulating Factor) in combination with one or more of:

[0019] collecting autologous cells and or plasma using a collection device and re-infusing the previously collected autologous cells and or blood plasma; and

[0020] re-infusing the previously collected autologous cells and or blood plasma together with allogeneic(healthy donor) cells and blood plasma.

[0021] Treating diseases of ageing and diseases of ageing with immunosenescence and immune dysfunction and inflammation by stem cell activation using a method of delivering precise magnetic field patterns which agree with the body's own natural magnetic field patterns, utilizing an instrument capable of routine clinical therapy use, in combination with healthy donor allogeneic cells and/or blood plasma.

[0022] Treating diseases of ageing and diseases of ageing with immunosenescence and immune dysfunction and inflammation by stem cell activation using a method of delivering precise magnetic field patterns which agree with the body's own natural magnetic field patterns, utilizing an instrument capable of routine clinical therapy use, in combination with stem cell mobilization agents, G-CSF, Granulocyte Colony Stimulating Factor)

[0023] Treating diseases of ageing and diseases of ageing with immunosenescence and

immune dysfunction and inflammation by stem cell activation using a method of delivering precise magnetic field patterns which agree with the body's own natural magnetic field patterns, utilizing an instrument capable of routine clinical therapy use.

[0024] Treating diseases of ageing and diseases of ageing with immunosenescence and immune dysfunction and inflammation by stem cell activation using a method of delivering precise magnetic field patterns which agree with the body's own natural magnetic field patterns, utilizing an instrument capable of routine clinical therapy use, in combination with autologous cells and/or blood plasma.

[0025] Other objects and advantages of this invention will become apparent from the following description taken in conjunction with any accompanying drawings wherein are set forth, by way of illustration and example, certain embodiments of this invention. Any drawings contained herein constitute a part of this specification and include exemplary embodiments of the present invention and illustrate various objects and features thereof.

BRIEF DESCRIPTION OF THE FIGURES

[0026] Figure 1 illustrates common electronic-generated signals;

[0027] Figure 2 illustrates the relationship between aging, immunosenescence, inflammation and disease states;

[0028] Figure 3 illustrates a SQUID installation at Vanderbilt University;

[0029] Figure 4 illustrates natural complex biologic waveforms found in the body;

[0030] Figure 5 illustrates natural complex biologic waveforms of a frog nerve;

[0031] Figure 6 illustrates the concept of extracting, analyzing and deriving waveforms;

[0032] Figure 7 illustrates magnetic field generation by current flow;

[0033] Figure 8 illustrates use of a solenoid design field generator;

- [0034] Figure 9 illustrates a toroidal magnetic field applicator;
- [0035] Figure 10 illustrates a Helmholtz coil field applicator;
- [0036] Figure 11 illustrates a planar field applicator;
- [0037] Figure 12 illustration of normal versus senescent cell morphology;
- [0038] Figure 13 illustrates an effect of the inventive method on regulation of levels of inflammatory and non-inflammatory markers in Patient 1;
- [0039] Figure 14 illustrates an effect of the inventive method on regulation of levels of inflammatory and non-inflammatory markers in Patient 2;
- [0040] Figure 15 illustrates an effect of the inventive method on regulation of levels of natural killer cells in Patient 2;
- [0041] Figure 16 illustrates an effect of the inventive method on regulation of levels of inflammatory and non-inflammatory markers in Patient 3;
- [0042] Figure 17 illustrates an effect of the inventive method on regulation of levels of naïve T cells in Patient 4;
- [0043] Figure 18 illustrates an effect of the inventive method on regulation of levels of Central Memory T cells and Natural Killer cell activity in Patient 4;
- [0044] Figure 19 illustrates an effect of the inventive method on regulation of levels of naïve T cells in Patient 5;
- [0045] Figure 20 illustrates an effect of the inventive method on regulation of levels of Central Memory T cells in Patient 5;
- [0046] Figure 21 illustrates an effect of the inventive method on regulation of levels of Natural Killer cell activity in Patient 5;
- [0047] Figure 22 illustrates an effect of the inventive method on regulation of levels of Natural Killer cell activity in Patient 6;

[0048] Figure 23 illustrates an effect of the inventive method on regulation of levels of Natural Killer cell activity in Patient 7;

[0049] Figure 24 is a SPECT Scan which demonstrates improvement in the patients Neurodegenerative Disease following treatment in accordance with the invention;

[0050] Figure 25 illustrates an effect of the inventive method on regulation of levels of Natural Killer cell activity in Patient 8;

[0051] Figure 26 illustrates an effect of the inventive method on regulation of levels of Total B cells in Patient 9;

[0052] Figure 27 illustrates an effect of the inventive method on regulation of levels of Memory T cells in Patient 9;

[0053] Figure 28 illustrates an effect of the inventive method on regulation of levels of naïve T cells in Patient 9;

[0054] Figure 29 illustrates an effect of the inventive method on regulation of levels of Natural Killer cell activity in Patient 9.

DETAILED DESCRIPTION OF THE INVENTION

[0055] Aging is associated with alterations of the immune system including impairments in innate immunity, T-lymphopoiesis and B-lymphopoiesis and these impairments contribute to immunosenescence in affected individuals.

[0056] An altered differentiation capacity of Hematopoietic stem cells (HSCs) has been causally linked to a reduction in lymphopoiesis during aging in mice and in man. Whole genome expression analyses indicated that HSC intrinsic alterations in gene expression contribute to this phenotype. The pool of HSCs comprises different HSC subpopulations that are biased toward myeloid or lymphoid differentiation. There is evidence that upon aging myeloid- biased HSCs are maintained, whereas lymphoid-biased HSCs get

lost. These result in the imbalance in myelolymphopoiesis occurring with aging. The molecular causes of this age-associated selection of HSC subpopulations remain to be delineated.

[0057] Accumulation of DNA damage has been associated with aging of HSCs in both mice and man. Moreover, studies on telomerase knockout mice (Terc ^{-/-}) revealed evidence that chronic DNA damage signaling in response to telomere dysfunction leads to an acceleration of hematopoietic skewing with a strong decrease in lymphopoiesis involving both cell intrinsic checkpoints and alterations in the blood circulatory environment. HSC aging contributes to impairments in early lymphoid lineage differentiation. This process associates with a selective increase of myeloid-competent HSCs and a decrease in lymphoid-competent HSCs during aging. This age-associated skewing in the maintenance of subpopulations of HSCs contributes to defects in lymphopoiesis and decreasing immune function during aging. Molecular mechanisms that can induce stem cell aging include the accumulation of DNA damage and telomere dysfunction and it is possible that stem cell intrinsic checkpoint as well as alteration in the stem cell environment (niche and systemic environment) can contribute to age-dependent selection of HSC subpopulations. The selective survival of distinct subpopulations of HSCs also contributes to the development of malignancies in the hematopoietic system and the selective maintenance of myeloid-competent HSCs enhances the risk of mutation accumulation in the myeloid lineage thereby leading to an increase of myelo-proliferative diseases during aging. The loss of lymphoid-competent HSCs may induce lymphoid lineage derived malignancies by impairing proliferative competition in lymphoid progenitor cell niches. Along these lines it has been shown that age-associated impairments in hematopoietic progenitor cell proliferation select for an outgrowth of malignant clones. In contrast to the possible influences on tumor promotion, it is conceivable that the depletion of HSC subpopulations could serve as a tumor

suppressor mechanism involved in the depletion of damaged HSCs. Cell surface marker combinations can subdivide human hematopoietic cells into different subpopulations which can also be subdivided into lymphoid-competent and/or myeloid-competent subpopulations during human aging. The stepwise process of the lymphoid differentiation of multipotent hematopoietic stem cells (HSCs) in human bone marrow has been assumed to begin with expression of the cell surface antigen CD10 (CALLA or MME) on CD34+ progenitors, based on the finding that CD 10+ progenitors lack myeloid and erythroid potential but are able to generate all lymphoid lineages. However, subsequent studies have shown that CD34+, CD10+ cells, even those without expression of lineage markers (Lin-: CD3-, CD14-, CD15-, CD19-, CD56-, CD235a-), show a strong bias toward B cell potential with relatively little T cell or natural killer (NK) cell potential. CD34+, Lin-, CD 10+ cells that lack expression of CD24 are precursors of the CD34+, Lin-CD 10+CD24+ population but nonetheless show molecular evidence of commitment to the B cell lineage, with expression of several B cell-specific genes. L-selectin (CD62L) is expressed on lymphocytes and mediates homing to peripheral lymphoid organs. Studies have reported that upregulation of CD62L expression on c-Kit+Lin-Sca-1+ mouse bone marrow cells correlates with loss of erythroid and megakaryocyte potential and efficient thymic engraftment. In the progenitor hierarchy of the lymphoid commitment of human cells, a stage of lymphoid priming that precedes commitment to the B lymphoid lineage, either before or independently of CD10 expression is a CD34+,Lin-, CD10- progenitor subpopulation in human bone marrow that has high expression of CD62L and that is devoid of clonogenic myeloid or erythroid potential. In stromal cultures, these cells are able to generate B cells, NK cells and T cells, as well as monocytic and dendritic cells.

[0058] Aging is associated with alterations of the immune system including impairments in T-lymphopoiesis and B-lymphopoiesis and these impairments contribute to

immunosenescence in affected individuals. Immunosenescence with immune dysfunction and increased inflammation is a primary cause of aging and diseases such as anemia of chronic diseases, autoimmune disorders, cancer, cardiovascular diseases, infection, metabolic diseases, neurodegenerative diseases, and protein energy malnutrition.

[0059] Natural magnetic field waveforms have been discovered associated with biologic processes ever since the discovery and development of the SQUID (Figure 3 "SQUID installation at Vanderbilt University"). These waveforms appear in complex patterns such as shown here in (Figure 4 "Natural complex biologic waveforms found in the body").

[0060] Although there is nothing new about the measurement and recording of these waveform the application of these waveforms in a useful clinic device has just recently been possible through the advance of modern electronics. Thus the method of identification, extraction and isolating and then delivering those magnetic field patterns in a therapeutically effective manner is a primary objective of the invention (Figure 6 "Extract, analyze and derive").

[0061] Researchers have theorized since the late 1960's that the information content of a magnetic field waveform is received and recognized by the body (if delivered in a specific manner) and is useful for therapeutic effects. The modern bone growth stimulator is one example of such a device, whereby this device has proven useful in medical applications to enhance the repair and growth of bone tissue. The method and device to deliver according to the method should, in theory, prove to be more effective in delivering medical therapy to a patient.

[0062] Natural biologic waveforms have been measured for many body processes. The graphic (Figure 4 "Natural complex biologic waveforms found in the body") describes several of these processes. More recently, sophisticated and sensitive recording technology has been used to record biologic processes with even greater sensitivity, such as the firing of a single nerve axon. Biologic waveforms have also been associated with specific diseases and inflammatory processes

which cause activation of stem cells local to the disease site and being mobilized from the bone marrow into the blood to reach the disease site. These waveforms are speculated to have an association with the body's natural healing processes. Researchers have also speculated that external electromagnetic signals applied to the body are ignored unless they are:

- Damaging signals, such as ionizing signals (e.g. X-Ray)
- Benign signals which affect a site of injury (e.g. diseases of ageing and diseases of ageing with immunosenescence and immune dysfunction and inflammation).

[0063] The device proposed herein is to deliver natural biologic waveform electromagnetic fields to a site of injury, and so function more effectively than the signals used previously in combination with methods of enhancing the concentration of stem cells at the disease sites. The method is the actual use of these natural biologic waveforms in the generation and delivery of those waveforms suited to a particular injury.

Process to capture, store and replicate biologic waveforms

[0064] The entire process begins with:

- Discovery of the biologic signal
- Isolation of the repair signal
- Storage of the repair signal
- Generation of the repair signal
- Delivery of the repair signal
- Conformity to a specified protocol

Discovery of the biologic signal

[0065] The discovery process begins with a known pathologic condition. For example, a cancer has well-understood biologic processes at work which serve to repair the immune dysfunction and immune senescence and inflammation. These processes all involve the generation and emitting of natural biologic waveforms.

[0066] The discovery, therefore, begins with a patient who has a known condition and a sensitive measurement device, known as a SQUID (Superconductive QUantum Interference Device), to detect and measure the condition waveforms. This device, or a representative of the device, is shown in Figure 3 "SQUID installation at Vanderbilt University". The waveforms generated by the body or biologic organisms have certain specific characteristics. Examples of certain waveforms are shown in Figure 5 "Natural complex biologic waveforms of a frog nerve".

[0067] The measurement of the natural biologic waveform caused by the underlying pathologic condition is facilitated by the SQUID apparatus, which is routinely used for measuring those types of waveforms. Conversely, the body may emit certain natural biologic waveforms that are associated with the normal biologic function. That is, those waveforms are captured from a healthy subject.

Isolation of the repair signal

[0068] The natural biologic waveform of the patient target pathology and the patient injury-free target are expected to differ in certain characteristics. In fact, Romanian researchers have reported in the literature that these signals do indeed exist and can be isolated. The isolation process may take place by digitizing those waveforms, analyzing and then performing certain digital operations on the patterns, using pattern-recognition or other graphical technology. Isolation of the waveforms is a straightforward procedure, by which measurements are taken of:

- A healthy subject
- A subject with diseases of ageing and diseases of ageing with immunosenescence and immune

dysfunction and inflammation

[0069] Each measurement is captured and digitized using mechanical or electrical conversion means and placed into a common file format. The procedure to further isolate the suspected natural biologic signal is a process whereby a comparison of the two waveforms yields a "difference" waveform (see "Extract, analyze and derive" on page 10), which is then presented as the suspected biologic waveform contributing to the healing process.

[0070] The original source waveforms, e.g. the "Disease" waveform, and The "Normal" waveform are used as reference waveforms, in a study to compare the relative effectiveness of those waveforms against the "difference" waveform.

Store of the repair signal

[0071] The final selected repair signal is then stored in electronic form, typically in a digitized fashion, or it may be stored in printed graphical form. This may use a common flat-file or relational database for the electronic storage medium.

Generation of the repair signal

[0072] The stored electrical signal pattern is then re-generated in a device (here called a Modulator) which then powers an external applicator. The re-generation of these electrical patterns may take place in a number of ways:

- Using an internal digital look-up table
- Using an internal derived equation, which is then solved
- Using a series of signals, such as is found in a

fourier series Once the basic repair signal is then re-generated, it is then modulated as to:

- Frequency
- Intensity

- Duty Cycle

by the generating device. This final signal is then amplified and prepared for delivery to the patient or subject.

Delivery of the repair signal

[0073] The repair signal, having been captured, stored, processed and modulated, is now ready to be delivered. Magnetic fields are delivered by passing a current through a wire loop assembly of various forms, as is generically illustrated in Figure 7. These types may be:

- A solenoid
- A toroid
- A planar (or flat) coil
- A Helmholtz coil

[0074] The solenoid magnetic field pattern may be seen in Figure 8 "Use of solenoid design". The toroid magnetic field patterns may be seen in Figure 9 "Toroidal magnetic field applicator". The planar coil design may be seen on Figure 11 "Planar field applicator". The Helmholtz design may be seen on Figure 10 "Helmholtz coil field applicator".

[0075] Prior to delivering the current to the various types of applicators, the Modulator power must amplify the stored signals to the level whereby they may energize the coils. This amplification may require the expenditure of power levels of up to 500 watts in the case of linear amplification, or may be significantly lower in the case of a switching- type (digital) design.

Conformity to a specified protocol

[0076] The conformity to a specified protocol may require specifying certain types of stage treatment procedures. These procedures typically require:

- A certain time of administered exposure
- A certain time between exposures
- A certain exposure dose
- A certain dose design
- A certain dose duty cycle

[0077] For example, the administered exposure may require stages exposures of 30, 60 or 90 minute exposures over the course of several days. The exposures may be staged according to a certain delay between each exposure. The dosage may be adjusted up or down according to the needs of the subject. The dose design itself is specified, according to a selection of one or several types of waveforms which are stored in the machine, and the dose may have to regulated according to a certain duty cycle.

[0078] Various embodiments of the invention are illustrated in the following examples.

Example 1

[0079] This protocol will assess the efficacy of stem cell activation with stem cell mobilization using a granulocyte colony stimulating factor (g-csf) and autologous stem cell rich plasma to improve the levels of anti-aging bio-markers in the recipients

[0080] .They will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly for 12 months. Assessments measuring the anti-aging bio-markers and clinical markers of aging as reported via a 13-item Clinical evaluation will be done at baseline before starting treatments and will be repeated every 3 months during

the treatment period, and every three months for 12 additional months after the completion of the treatment period. Half of the recipients will be randomly selected to also receive Autologous Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months in the second 12 month study period. The Autologous plasma donors are also plasma recipients. Autologous Donors will be dosed with G-CSF for three days prior to undergoing plasmapheresis. This will stimulate a significantly increased number of stem cells in the plasma.

The yield of stem cell rich plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots. Samples will be frozen until such time as the same Plasma Recipient is being treated.

[0081] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[0082] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Safety Assessments:

[0083] Baseline physical exam, blood chemistry will be done on all recipients at baseline. These will be repeated after 3, 6, 9, and 12 months of treatment, and quarterly for 1 year following the completion of the last treatment.

[0084] Safety and tolerability will be monitored through continuous reporting of adverse events.

Example 2

[0085] This protocol will assess the efficacy of stem cell activation with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) and autologous stem cell rich plasma to improve the levels of anti-aging bio-markers in the recipients

[0086] Patients will be treated with G-CSF stem cell mobilization and Autologous Stem Cell Rich Plasma to improve the levels of anti-aging biomarkers in the recipients.

Primary Objective:

[0087] During treatment, patients will continue to be evaluated for improvement in the levels of anti-aging biomarkers in the recipients during treatment for one year and for one year following the end of treatment.

Secondary Objectives:

[0088] Improvement of clinical markers of aging as reported via a 13-item Clinical evaluation

Study Design:

[0089] The study is classified as Phase I/II since it will be using a long-standing safe 'drug' (G-CSF; Autologous plasma; ergo Phase II) but in a novel way (ergo Phase I).

[0090] They will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly for 12 months. On the first of the 5-7 day period they will also receive Autologous Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[0091] Autologous Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. The yield of stem cell rich plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots. Samples will be frozen until such time the same Plasma Recipient is being treated.

[0092] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

- 1) Immunosenescence Panel
- 2) CCL 11
- 3) TGFBeta 1 growth factor
- 4) Nuclear Factor kappa beta (NFkB)
- 5) DHEA-S

- 6) Plasma Insulin
- 7) Telomere length
- 8) Cytokine Multiplex 18

[0093] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Example 3

[0094] This protocol will assess the efficacy of stem cell activation with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) and autologous stem cell rich plasma in combination with infusing stem-cell rich plasma from ABO-matched cord blood healthy allogeneic donors to improve the levels of anti-aging bio-markers in the recipients.

[0095] Patients will be treated with G-CSF stem cell mobilization and Autologous Stem Cell Rich Plasma in combination with Infusing Stem-Cell Rich Plasma from ABO-matched Healthy Cord Blood Allogeneic donors to improve the levels of anti-aging biomarkers in the recipients.

[0096] Adults who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of G-CSF stem cell mobilization as well as Autologous Stem Cell Rich Plasma and Infusions of Stem-Cell Rich Plasma from ABO-matched Healthy Cord Blood Allogeneic Donors will be treated monthly for a one-year period to determine if there is any improvement in those markers.

[0097] They will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly for 12 months. On the first of the 5-7 day period they will also receive Autologous Stem Cell Rich plasma in aliquots of 50 ml monthly for 12

months. After the end of the daily G-CSF for the 5-7 day period they will also receive Allogeneic Cord Blood Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[0098] Autologous Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasmapheresis. This will stimulate a significantly increased number of stem cells in the plasma.

[0099] The yield of stem cell rich plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots. Samples will be frozen until such time as the same Plasma Recipient is being treated.

Allogeneic Cord Blood Plasma Donors:

[00100] All donations of cord blood plasma will be taken from the plasma extracted and normally discarded from healthy babies cord blood collected from the delivered placenta at the time of birth and whose stem cells are to be cryopreserved for the baby's personal use at a later date.

[00101] The yield of plasma expected from the cord blood is about 55 ml. The plasma will be divided into aliquots of 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time as a Plasma Recipient of the same ABO/Rh type is being treated.

[00102] Based on the average prevalence of ABO types, it is anticipated that the number of donors needed will be approximately 12 to 15 per recipient.

[00103] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00104] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Example 4

[00105] Recipients will be treated with G-CSF stem cell mobilization and Autologous Stem Cell Rich Plasma in combination with Infusing Stem-Cell Rich Plasma from ABO-matched Healthy Allogeneic Donors to improve the levels of anti-aging biomarkers in the recipients.

[00106] **Recipients** who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of G-CSF stem cell mobilization as well as Autologous Stem Cell Rich Plasma and Infusions of Stem-Cell Rich Plasma from ABO-

matched Healthy Allogeneic Donors will be treated monthly for a one-year period to determine if there is any improvement in those markers. I).

[00107] They will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days .monthly for 12 months. On the first of the 5-7 day period they will also receive Autologous Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. After the end of the daily G-CSF for the 5-7 day period they will also receive Allogeneic Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00108] Autologous Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. The yield of stem cell rich plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots. Samples will be frozen until such time the same Plasma Recipient is being treated. .

Allo2eneic Stem Cell Rich Plasma Donors :

[00109] The other group making up the study population will be Allogeneic plasma donors. Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. Donors will have Infectious Disease testing repeated on the day of pheresis.

[00110] The yield of plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time as a Plasma Recipient of the same ABO/Rh type is being treated. Donors will be restricted to males and nulliparous females to avoid the presence of cytotoxic lymphocyte and granulocyte antibodies.

[00111] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00112] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Example 5

[00113] This protocol will assess the efficacy of stem cell activation with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) and in combination with

infusing stem-cell rich plasma from ABO-matched healthy allogeneic donors to improve the levels of anti-aging bio-markers in the recipients

[00114] Patients will be treated with G-CSF stem cell mobilization and in combination with Infusing Stem-Cell Rich Plasma from ABO-matched Healthy Allogeneic Donors will be monitored for the improvement of levels of anti-aging biomarkers in the recipients.

[00115] **Recipients** who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of G-CSF stem cell mobilization and Infusions of Stem-Cell Rich Plasma from ABO-matched Healthy Allogeneic Donors will be treated monthly for a one-year period to determine if there is any improvement in those markers. Recipients will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly for 12 months. After the end of the daily G-CSF for the 5-7 day period they will also receive Allogeneic Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00116] **Allogeneic Stem Cell Rich Plasma Donors:**

Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. The yield of plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until

such time as a Plasma Recipient of the same ABO/Rh type is being treated. Donors will be restricted to males and nulliparous females to avoid the presence of cytotoxic lymphocyte and granulocyte antibodies.

[00117] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00118] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Example 6

[00119] This protocol will assess the efficacy of stem cell activation with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) and autologous stem cell rich plasma in combination with infusing stem-cell rich plasma from ABO-matched healthy allogeneic donors to improve the levels of anti-aging bio-markers in the recipients.

[00120] Patients will be treated with G-CSF stem cell mobilization and Autologous Stem Cell Rich Plasma in combination with Infusing Stem-Cell Rich Plasma from ABO-matched Healthy Allogeneic Donors in the hope of improving the levels of anti-aging biomarkers in the recipients.

[00121] Recipients who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of G-CSF stem cell mobilization as well as Autologous Stem Cell Rich Plasma and Infusions of Stem-Cell Rich Plasma from ABO-matched Healthy Allogeneic Donors will be treated monthly for a one-year period to determine if there is any improvement in those markers.

Recipients will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly for 12 months. On the first of the 5-7 day period they will also receive Autologous Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. After the end of the daily G-CSF for the 5-7 day period they will also receive Allogeneic Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00122] Autologous Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. The yield of stem cell rich plasma expected from the pheresis is about 2 liters. The plasma will be

divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots. Samples will be frozen until such time as the same Plasma Recipient is being treated.

Allogeneic Stem Cell Rich Plasma Donors:

Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. The yield of plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time as a Plasma Recipient of the same ABO/Rh type is being treated. Donors will be restricted to males and nulliparous females to avoid the presence of cytotoxic lymphocyte and granulocyte antibodies.

[00123] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00124] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Example 7

[00125] This protocol will assess the efficacy of stem cell activation with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) in combination with infusing stem-cell rich plasma from ABO-matched cord blood allogeneic donors to improve the levels of anti-aging bio-markers in the recipients.

[00126] Patients will be treated with G-CSF stem cell mobilization in combination with Infusing Stem-Cell Rich Plasma from ABO-matched Cord Blood Allogeneic Donors will improve the levels of anti-aging biomarkers in the recipients.

[00127]

[00128] **Recipients** who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of G-CSF stem cell mobilization as well as Infusions of Stem-Cell Rich Plasma from ABO-matched Cord Blood Allogeneic Donors will be treated monthly for a one-year period to determine if there is any improvement in those markers.

[00129] Recipients will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly. After the end of the daily G-CSF for the 5-7 day period, they will receive Allogeneic Cord Blood Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

Cord Blood Plasma Donors:

[00130] All donations of cord blood plasma will be taken from the plasma extracted and normally discarded from healthy babies cord blood collected from the delivered placenta at the time of birth and whose stem cells are to be cryopreserved for the baby's personal use at a later date. The yield of plasma expected from the cord blood is about 55 ml. The plasma will be divided into aliquots of about 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time as a Plasma Recipient of the same ABO/Rh type is being treated.

[00131] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFβ1 growth factor
4. Nuclear Factor kappa beta (NFκB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00132] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Example 8

[00133] This protocol will assess the efficacy of stem cell activation using the recipients own natural magnetic field patterns in combination with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) and autologous stem cell rich plasma to improve the levels of anti-aging bio-markers in the recipients.

[00134] Patients will be treated with G-CSF stem cell mobilization and Autologous Stem Cell Rich Plasma in Combination with administration of precise natural Magnetic fields to improve the levels of anti-aging biomarkers in the recipients and improvement of clinical markers of aging as reported via a 13-item Clinical evaluation

[00135] **Recipients** who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of G-CSF stem cell mobilization as well as Autologous Stem Cell Rich Plasma and administration of their own Precise Natural Magnetic fields will be treated monthly for a one-year period to determine if there is any improvement in those markers.

[00136] Recipients will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly for 12 months. On the first of the 5-7 day period they will also receive Autologous Stem Cell Rich plasma in aliquots of 50 ml and precise magnetic field patterns which will agree with their body's own natural magnetic field patterns monthly for 12 months. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00137] **Autologous** Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma.

The yield of stem cell rich plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots.

Samples will be frozen until such time as the same Plasma Recipient is being treated.

[00138] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00139] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Example 9

[00140] This protocol will assess the efficacy of stem cell activation using the recipients own natural magnetic field patterns in combination with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) and autologous stem cell rich plasma in combination with infusing stem-cell rich plasma from ABO-matched healthy allogeneic donors to improve the levels of anti-aging bio-markers in the recipients.

[00141] Patients will be treated with G-CSF stem cell mobilization and Autologous Stem Cell Rich Plasma in combination with Infusing Stem-Cell Rich Plasma from ABO-matched Healthy Allogeneic Donors in Combination with administration of Precise natural Magnetic fields to improve the levels of anti-aging biomarkers in the recipients.

[00142] **Recipients** who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of G-CSF stem cell mobilization as well as Autologous Stem Cell Rich Plasma and Infusions of Stem-Cell Rich Plasma from ABO-matched Healthy Allogeneic Donors and administration of their own Precise Natural Magnetic fields monthly for a one-year period will be treated to determine if there is any improvement in those markers. The recipients will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly for 12 months. On the first day of the 5-7 day period they will also receive Autologous Stem Cell Rich plasma in aliquots of 50 ml and precise magnetic field patterns which will agree with their body's own natural magnetic field patterns monthly for 12 months. After the end of the daily G-CSF for the 5-7 day period they will also receive Allogeneic Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00143] Autologous Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasmapheresis. This will stimulate a significantly increased number of stem cells in the plasma. The yield of stem cell rich plasma expected from the pheresis is about 2 liters. The plasma will be divided into

aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots. Samples will be frozen until such time the same Plasma Recipient is being treated.

Allo2eneic Stem Cell Rich Plasma Donors :

Allogeneic plasma donors will be healthy, young adults (age 30 or less) with no major medical diagnoses. A physical exam and standard blood chemistry and CBC will determine their eligibility. An ABO/Rh typing, hemoglobinopathy testing and antibody panel will be done on each donor, as well as Infectious Disease testing. Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. Donors will have Infectious Disease testing repeated on the day of pheresis. The yield of plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time a Plasma Recipient of the same ABO/Rh type is being treated. Donors will be restricted to males and nulliparous females to avoid the presence of cytotoxic lymphocyte and granulocyte antibodies.

[00144] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFβ1 growth factor
4. Nuclear Factor kappa beta (NFκB)
5. DHEA-S
6. Plasma Insulin

7. Telomere length

8. Cytokine Multiplex 18

[00145] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation. Baseline physical exam, blood chemistry will be done on all recipients at baseline. These will be repeated after 3, 6, 9, and 12 months of treatment, and quarterly for 1 year following the completion of treatment.

Example 10

[00146] This protocol will assess the efficacy of stem cell activation using natural magnetic field patterns in combination with infusing stem-cell rich plasma from ABO-matched healthy allogeneic donors to improve the levels of anti-aging bio-markers in the recipients.

[00147] Recipients will be treated with plasma from healthy allogeneic donors who have had stem cell mobilization in combination with administration of Precise natural Magnetic fields to improve the levels of anti-aging biomarkers in the recipients of the plasma.

[00148] Adults who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of donor plasma and on the same day they will also receive precise magnetic field patterns which will agree with their body's own natural magnetic field patterns will be treated monthly for a one-year period to determine if there is any improvement in those markers. The recipients will receive ABO/Rh cross-matched plasma in aliquots of 50 ml monthly for 12 months. On the same day they will receive administration of Precise natural Magnetic fields. Assessments measuring the biomarkers will be done at

baseline and every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers(group of 13) will be collected at those same time points.

Plasma Donors:

[00149] All donors will be healthy, young adults (age 30 or less) with no major medical diagnoses. A physical exam and standard blood chemistry and CBC will determine their eligibility. An ABO/Rh typing and antibody panel will be done on each donor, as well as Infectious Disease testing. Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. Donors will have Infectious Disease testing repeated on the day of pheresis.

The yield of plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time as a Plasma Recipient of the same ABO/Rh type is being treated. Donors will be restricted to males and nulliparous females to avoid the presence of cytotoxic lymphocyte and granulocyte antibodies.

[00150] Recipients can be of any age or gender as long as they have no significant medical issues that would be contraindicated by the infusion of donor plasma and administration of precise magnetic field patterns which will agree with their body's own natural magnetic field patterns. This will be determined by a physical exam and standard blood chemistry and CBC. Their ABO/Rh type will also be tested.

[00151] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00152] There will also be clinical markers (group of 13) evaluated at baseline and at quarterly visits during physical exam that will provide secondary efficacy data for evaluation. Baseline physical exam, blood chemistry and infectious disease markers will be done on all plasma donors and recipients at baseline. These apart from infectious disease markers will be repeated after 3, 6, 9, and 12 months of treatment, and quarterly for 1 year following the completion of treatment.

Example 11

[00153] This protocol will assess the efficacy of stem cell activation using natural magnetic field patterns in combination with infusing stem-cell rich plasma from ABO-matched healthy cord blood allogeneic donors to improve the levels of anti-aging biomarkers in the recipients.

[00154] Recipients will be treated with plasma from healthy Cord Blood allogeneic donors in combination with administration of Precise natural Magnetic fields to improve the levels of anti-aging biomarkers in the recipients of the plasma.

[00155] Recipients who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of Cord Blood donor plasma and on the same day they will also receive precise magnetic field patterns which will agree with their body's own natural magnetic field patterns will be treated monthly for a one-year period to determine if there is any improvement in those markers.

[00156] The recipients will receive ABO/Rh cross-matched Cord Blood plasma in aliquots of 50 ml monthly for 12 months. On the same day they will receive administration of Precise natural Magnetic fields. Assessments measuring the biomarkers will be done at baseline and every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00157] All donations of cord blood plasma will be taken from the plasma extracted and normally discarded from healthy babies cord blood collected from the delivered placenta at the time of birth and whose stem cells are to be cryopreserved for the baby's personal use at a later date. The baby and mother will have no major medical diagnoses and the mother will sign an informed consent form and agree to donation of cord blood plasma. A physical exam and standard blood chemistry and CBC will determine their eligibility. An ABO/Rh typing and antibody panel will be done on each donor, as well as Infectious Disease testing.

[00158] The yield of plasma expected from the cord blood is about 55 ml. The plasma will be divided into aliquots of 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time as a Plasma Recipient of the same ABO/Rh type is being treated. Based on the average prevalence of ABO types, it is anticipated that the number of donors needed will be approximately 12 to 15 per recipient.

[00159] Recipients can be of any age or gender as long as they have no significant medical issues that would be contraindicated by the infusion of donor plasma and administration of precise magnetic field patterns which will agree with their body's own natural magnetic field patterns. This will be determined by a physical exam and standard blood chemistry and CBC. Their ABO/Rh type will also be tested.

[00160] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

There will also be clinical markers (group of 13) evaluated at baseline and at quarterly visits during physical exam that will provide secondary efficacy data for evaluation. Baseline physical exam, blood chemistry and infectious disease markers will be done on all plasma donors and recipients at baseline. These apart from infectious disease markers will be repeated after 3, 6, 9, and 12 months of treatment, and quarterly for 1 year following the completion of treatment.

Example 12

[00161] This protocol will assess the efficacy of stem cell activation using the recipients own natural magnetic field patterns in combination with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) to improve the levels of anti-aging bio-markers in the recipients.

[00162] Patients will be treated with G-CSF stem cell mobilization in combination with administration of Precise natural Magnetic fields to improve the levels of anti-aging biomarkers in the recipients.

[00163] Recipients who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of G-CSF stem cell mobilization and administration of their own Precise natural Magnetic fields ill be treated monthly for a one-year period to determine if there is any improvement in those markers.

The recipients will receive daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for 5-7 days monthly for 12 months. On the same day they will also receive precise magnetic field patterns which will agree with their body's own natural magnetic field patterns. Assessments measuring the biomarkers will be done at baseline before starting treatments and will be repeated every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points._All recipients will have no major medical diagnoses which exclude them from undergoing G-CSF stem cell mobilization and administration of Precise natural Magnetic fields which agree with the recipients own natural magnetic field patterns utilizing an instrument capable of routine clinical therapy use and assessment of response. . A physical exam, standard blood chemistry ,CBC and hemoglobinopathy testing will be done.

[00164] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same

evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFβ1 growth factor
4. Nuclear Factor kappa beta (NFκB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00165] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation. Baseline physical exam, blood chemistry will be done on all recipients at baseline. These will be repeated after 3, 6, 9, and 12 months of treatment, and quarterly for 1 year following the completion of treatment.

Example 13

[00166] This protocol will assess the efficacy of stem cell activation using the recipients own natural magnetic field patterns to improve the levels of anti-aging biomarkers in the recipients.

[00167] Recipients will be treated by administration of Precise natural Magnetic fields to improve the levels of anti-aging biomarkers.

[00168] Evaluation of the improvement in the levels of anti-aging biomarkers in the recipients during treatment for one year and for one year following the end of treatments.

Recipients who are interested in assessing their level of anti-aging biomarkers, and who are

willing to receive stem cell activation using Precise natural Magnetic fields will be treated monthly for a one-year period to determine if there is any improvement in those markers. Recipients will receive precise magnetic field patterns which will agree with their body's own natural magnetic field patterns monthly for 12 months. Assessments measuring the biomarkers will be done prior to the start of treatment and thereafter every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

All recipients will have no major medical diagnoses which exclude them from receiving Precise natural Magnetic fields which agree with the recipients own natural magnetic field patterns utilizing an instrument capable of routine clinical therapy use and assessment of response. A physical exam and standard blood chemistry ,CBC will determine their eligibility.

[00169] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00170] There will also be clinical markers (group of 13) evaluated at baseline and at quarterly visits during physical exam that will provide efficacy data for evaluation. Baseline physical exam, blood chemistry and biomarkers markers will be done on all recipients at baseline. These will be repeated after 3, 6, 9, and 12 months of treatment, and quarterly for 1 year following the completion of treatment.

Example 14

[00171] This protocol will assess the efficacy of stem cell activation using the recipients own natural magnetic field patterns in combination with infusing autologous stem-cell rich plasma to improve the levels of anti-aging bio-markers in the recipients.

[00172] Recipients will be treated with autologous stem cell rich plasma from recipients who have had stem cell mobilization in combination with administration of Precise natural Magnetic fields to improve the levels of anti-aging biomarkers in the recipients. Evaluation of the improvement in the levels of anti-aging biomarkers in the recipients during treatment for one year and for one year following the end of treatment.

Recipients who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of autologous stem cell rich donor plasma will be treated monthly for a one-year period to determine if there is any improvement in those markers.

[00173] Recipients will receive Autologous Stem Cell Rich plasma in aliquots of 50 ml monthly for 12 months. On the same day they will also receive precise magnetic field patterns which will agree with their body's own natural magnetic field patterns. Assessments measuring the biomarkers will be done every 3 months during the treatment period, and every

three months for 12 additional months after the completion of the treatment period.

Additionally, clinical markers will be collected at those same time points.:

[00174] All autologous donors who are also recipients will have no major medical diagnoses which exclude them from undergoing stem cell mobilization . A physical exam and standard blood chemistry, CBC ,and Infectious Disease marker testing will determine their eligibility. Their ABO/Rh type will also be tested.

[00175] Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. Donors will have Infectious Disease testing repeated on the day of pheresis. ABO/Rh type will also be tested.

[00176] The yield of stem cell rich plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots. Samples will be frozen until such time as the same Plasma Recipient is being treated.

[00177] Administration of Precise natural Magnetic fields which agree with the recipients own natural magnetic field patterns utilizing an instrument capable of routine clinical therapy use and response assessment

[00178] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)

5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00179] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Baseline physical exam, blood chemistry and infectious disease markers will be done on all plasma donors and recipients at baseline. These apart from infectious disease markers will be repeated after 3, 6, 9, and 12 months of treatment, and quarterly for 1 year following the completion of treatment.

Example 15

[00180] This protocol will assess the efficacy of infusing stem-cell rich plasma from ABO-matched healthy donors to improve the levels of anti-aging bio-markers in the recipients. Recipients will be treated with plasma from healthy donors who have had stem cell mobilization to improve the levels of anti-aging biomarkers in the recipients of the plasma.

[00181] Recipients who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of donor plasma will be treated monthly for a one-year period to determine if there is any improvement in those markers.

[00182] The recipients will receive ABO/Rh cross-matched plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done every 3 months during the treatment period, and every three months for 12 additional months after

the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00183] All donors will be healthy, young adults (age 30 or less) with no major medical diagnoses. A physical exam and standard blood chemistry and CBC will determine their eligibility. An ABO/Rh typing and antibody panel will be done on each donor, as well as Infectious Disease testing.

[00184] Donors will be dosed with G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously for three days prior to undergoing plasma pheresis. This will stimulate a significantly increased number of stem cells in the plasma. Donors will have Infectious Disease testing repeated on the day of pheresis.

[00185] The yield of plasma expected from the pheresis is about 2 liters. The plasma will be divided into aliquots of about thirty-six 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time as a Plasma Recipient of the same ABO/Rh type is being treated. Donors will be restricted to males and nulliparous females to avoid the presence of cytotoxic lymphocyte and granulocyte antibodies.

[00186] Recipients can be of any age or gender as long as they have no significant medical issues that would be contraindicated by the infusion of donor plasma. This will be determined by a physical exam and standard blood chemistry and CBC. Their ABO/Rh type will also be tested.

[00187] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11

3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00188] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Baseline physical exam, blood chemistry and infectious disease markers will be done on all plasma donors and recipients at baseline. These apart from infectious disease markers will be repeated after 3, 6, 9, and 12 months of treatment, and quarterly for 1 year following the completion of treatment.

Example 16

[00189] This protocol will assess the efficacy of infusing stem-cell rich plasma from A β O-matched cord blood donors to improve the levels of anti-aging bio-markers in the recipients.

[00190] Recipients will be treated with plasma from Cord Blood donors to improve the levels of anti-aging biomarkers in the recipients of the plasma.

Evaluation of the improvement in the levels of anti-aging biomarkers in the recipients during treatment for one year and for one year following the end of treatment.

Improvement of clinical markers of aging as reported via a 13-item Clinical evaluation.

[00191] **Recipients** who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive doses of donor plasma monthly for a one-year period will be treated to determine if there is any improvement in those markers.

[00192] The recipients will receive ABO/Rh cross-matched Cord Blood Plasma in aliquots of 50 ml monthly for 12 months. Assessments measuring the biomarkers will be done every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00193] All donations of cord blood plasma will be taken from the plasma extracted and normally discarded from healthy babies cord blood collected from the delivered placenta at the time of birth and whose stem cells are to be cryopreserved for the baby's personal use at a later date. The baby and mother will have no major medical diagnoses and the mother will sign an informed consent form and agree to donate plasma. A physical exam and standard blood chemistry and CBC will determine their eligibility. An ABO/Rh typing and antibody panel will be done on each donor, as well as Infectious Disease testing. The yield of plasma expected from the cord blood is about 55 ml. The plasma will be divided into aliquots of about 50 ml doses and 5 ml testing aliquots (for cross matching). Samples will be frozen until such time as a Plasma Recipient of the same ABO/Rh type is being treated.

[00194] Recipients can be of any age or gender as long as they have no significant medical issues that would be contraindicated by the infusion of donor plasma. This will be determined by a physical exam and standard blood chemistry and CBC. Their ABO/Rh type will also be tested.

[00195] There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same

evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFBeta 1 growth factor
4. Nuclear Factor kappa beta (NFkB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00196] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

Example 18

[00197] This protocol will assess the efficacy of G-CSF (granulocyte colony stimulating factor) mobilization of CD34+ peripheral blood stem cells to improve the levels of anti-aging bio-markers in the recipients.

[00198] Patients will be treated with stem cell mobilization factor to improve the levels of anti-aging biomarkers in the recipients. Evaluation of the improvement in the levels of anti-aging biomarkers in the recipients during treatment for one year and for one year following the end of treatment. Improvement of clinical markers of aging as reported via a 13-item Clinical evaluation.

[00199] **Recipients** who are interested in assessing their level of anti-aging biomarkers, and who are willing to receive G-CSF Mobilization for a one-year period will be treated to determine if there is any improvement in those markers.

[00200] The recipients will receive initially 3 cycles of G-CSF mobilization followed by 1 cycle at 3 months, 6 months, 9 months and 12 months. Each Cycle consists of daily G-CSF at a dose of 5 to 15ug/kg intravenously or subcutaneously per day administered subcutaneously for 3 to 7 days followed by 7 days off G-CSF with evaluation on 3 of the 7 off days. Assessments measuring the biomarkers will be done every 3 months during the treatment period, and every three months for 12 additional months after the completion of the treatment period. Additionally, clinical markers will be collected at those same time points.

[00201] Recipients who sign an informed consent form can be of any age or gender as long as they have no significant medical issues that would be contraindicated by the subcutaneous administration of G-CSF. This will be determined by a physical exam and standard blood chemistry and CBC. Their hemoglobinopathy screen will also be tested. There will be eight primary biomarker measurements evaluated at baseline, and after treatment for 3 months, 6 months, 9 months, and 12 months. These same evaluations will continue every three months for an additional year to examine the long-term effect of the treatment. These assessments are:

1. Immunosenescence Panel
2. CCL 11
3. TGFβ1 growth factor
4. Nuclear Factor kappa beta (NFκB)
5. DHEA-S
6. Plasma Insulin
7. Telomere length
8. Cytokine Multiplex 18

[00202] There will also be clinical markers (group of 13) evaluated at the quarterly visits during physical exam that will provide secondary efficacy data for evaluation.

[00203] In accordance with the protocol set forth in Example 18, several patients were treated. These patients, their conditions, and treatment will now be summarized with reference to Figures 13-29.

[00204] Regarding Patient 1, reference is made to Figure 13.

[00205] Diagnosis: Anemia, Chronic Disease(Chronic Obstructive Pulmonary Disease), Cardiovascular Disease (Chronic Heart Failure), Protein energy malnutrition, and Frailty.

[00206] In accordance with the basic protocol outlined in Example 18, a patient MR was treated with stem cell mobilization factor, G-CSF to improve the levels of anti-aging biomarkers in the patient. Evaluation of the levels of anti-inflammatory(IL-10) and inflammatory(TNF-alpha) aging biomarkers during treatment was carried out. An improvement in the level of the clinical marker of aging, IL-10 was reported. After 3 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of IL-10 improved from 4.54 to 9.489.

[00207] Evaluation of the improvement in the levels of inflammatory aging biomarkers during treatment was carried out. A decrease in the level of the inflammatory clinical marker of aging, TNF-alpha was reported. After 3 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of TNF-alpha went from 20.636 to 9.997.

[00208] While not wishing to be bound to any particular theory of operation, reduction in TNF-alpha and increase in IL-10 resulted in improvement in clinical symptoms of inflammation and immunosenescence associated with Anemia, Chronic Disease(Chronic Obstructive Pulmonary Disease,), Cardiovascular Disease (Chronic Heart Failure), Protein energy malnutrition, Frailty.

[00209] Regarding Patient 2, reference is made to Figures 14 and 15.

[00210] Diagnosis: Cancer, Protein energy malnutrition, Frailty.

[00211] In accordance with the basic protocol outlined in Example 18, a patient DP was treated with stem cell mobilization factor, G-CSF to improve the levels of aging biomarkers in the patient. Evaluation of the levels of anti-inflammatory(IL-10) and inflammatory(TNF-alpha) aging biomarkers during treatment was carried out. An improvement in the level of the clinical marker of aging, IL-10 was reported. After 3 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of IL-10 improved from 1.787 to 5.774.

[00212] Evaluation of the improvement in the levels of inflammatory aging biomarkers during treatment was carried out. A decrease in the level of the inflammatory clinical marker of aging, TNF-alpha was reported. After 3 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of TNF-alpha decreased from 11.072 to 7.243. After 3 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of Natural Killer Cells increased from 82 to 183.

[00213] While not wishing to be bound to any particular theory of operation, reduction in TNF-alpha and increase in IL-10 and Natural Killer Cells resulted in improvement in clinical symptoms of inflammation and immunosenescence associated with Cancer, Protein energy malnutrition, Frailty.

[00214] Regarding Patient 3, reference is made to Figure 16.

[00215] Diagnosis: Chronic disease, Neurodegenerative Disease, Frailty.

[00216] In accordance with the basic protocol outlined in Example 18, a patient JB was treated with stem cell mobilization factor, G-CSF to improve the levels of aging biomarkers in the patient. Evaluation of the levels of anti-inflammatory(IL-10) and inflammatory(TNF-alpha) aging biomarkers during treatment was carried out. An improvement in the level of the clinical marker of aging, IL-10 was reported. After 3 months of treatment, during which

time the patient was administered G-CSF as in Example 18, levels of IL-10 improved from 4.326 to 6.264.

[00217] Evaluation of the improvement in the levels of inflammatory aging biomarkers during treatment was carried out. A decrease in the level of the inflammatory clinical marker of aging, TNF-alpha was reported. After 3 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of TNF-alpha decreased from 9.469 to 3.832.

[00218] While not wishing to be bound to any particular theory of operation, reduction in TNF-alpha and increase in IL-10 resulted in improvement in clinical symptoms of inflammation and immunosenescence associated with Chronic disease, Neurodegenerative Disease ,Frailty.

[00219] Regarding Patient 4, reference is made to Figures 17 and 18.

[00220] Diagnosis: Chronic Metabolic Disease, Diabetes Mellitus Type 2, Frailty.

[00221] In accordance with the basic protocol outlined in Example 18, a patient RA was treated with stem cell mobilization factor, G-CSF to improve the levels of aging biomarkers in the patient. Evaluation of the improvement in the levels of aging, immunosenescence, immune dysfunction, and early lymphoid lineage differentiation biomarkers during treatment was carried out. An improvement with an increase in the level of the Naive CD4 and CD8 levels and a decrease in the Memory CD4 T cells and improvement in Natural Killer Cell Activity were reported. After 2 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of Naive CD4 improved from 4.88 to 12.46. Levels of Naive CD8 improved from 14.62 to 24.17.

[00222] Evaluation of the improvement was observed with a decrease in the levels of Memory CD4 T Cells from 80.37 to 68.52.

[00223] Also an improvement in Natural Killer Cell Activity was reported After 2 months of treatment, during which time the patient was administered G-CSF as in Example 18, wherein levels of Natural Killer Cell Activity CD4 improved from 9.46 to 18.36.

[00224] While not wishing to be bound to any particular theory of operation, reduction in Central Memory T cells and increase in Naive T cells and Natural Killer Cell Activity resulted in improvement in clinical symptoms of immune dysfunction, immunosenescence and impairment of early lymphoid differentiation associated with Chronic Metabolic Disease, Diabetes Mellitus Type 2, and Frailty.

[00225] Regarding Patient 5, reference is made to Figures 19-21.

[00226] Diagnosis: Chronic disease, Cancer (Waldenstroms Macroglobulinemia), Neurodegenerative Disease (Peripheral Neuropathy), Frailty.

[00227] In accordance with the basic protocol outlined in Example 18, a patient DS was treated with stem cell mobilization factor, G-CSF to improve the levels of aging biomarkers in the patient. Evaluation of the improvement in the levels of aging, immunosenescence, immune dysfunction, and early lymphoid lineage differentiation biomarkers during treatment was carried out. An improvement with an increase in the level of the Naive CD4 levels and a decrease in the Memory CD4 T cells and improvement in Natural Killer Cell Activity was reported. After 2 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of Naive CD4 improved from 12.08 to 42.47.

[00228] Evaluation of the improvement in immune dysfunction was further observed with a decrease in the levels of Memory CD4 T Cells decreased from 78.09 to 40.90.

[00229] Also an improvement in Natural Killer Cell Activity was reported. After 2 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of Natural Killer Cell Activity CD4 improved from 1.93 to 7.22.

[00230] While not wishing to be bound to any particular theory of operation, reduction in Central Memory T cells and increase in Naive T cells and Natural Killer Cell Activity resulted in improvement in clinical symptoms of immune dysfunction, immunosenescence and impairment of early lymphoid differentiation associated with Chronic disease, Cancer (Waldenstroms Macroglobulinemia), Neurodegenerative Disease (Peripheral Neuropathy), and Frailty.

[00231] Regarding Patient 6, reference is made to Figure 22.

[00232] Diagnosis: Chronic disease, Chronic Infection (Lyme's Disease) , Neurodegenerative Disease.

[00233] In accordance with the basic protocol outlined in Example 18, a patient SH was treated with stem cell mobilization factor, G-CSF to improve the levels of aging biomarkers in the patient. Evaluation of the improvement in the levels of aging, immunosenescence, and immune dysfunction biomarkers during treatment was carried out. An improvement with an increase in the level of Natural Killer Cells was reported. After 2 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of Natural Killer Cells improved from 268 to 623.

[00234] While not wishing to be bound to any particular theory of operation increase in Natural Killer Cells resulted in improvement in clinical symptoms of immune dysfunction, and immunosenescence associated with Chronic disease, Chronic Infection (Lyme's Disease) , and Neurodegenerative Disease.

[00235] Regarding Patient 7, reference is made to Figures 23 and 24.

[00236] Diagnosis: Chronic disease, Chronic Fatigue Syndrome, Neurodegenerative Disease, Autoimmune Disease(Type 1Diabetes Mellitus).

[00237] In accordance with the basic protocol outlined in Example 18, a patient LS was treated with stem cell mobilization factor, G-CSF to improve the levels of aging

biomarkers in the patient. Evaluation of the improvement in the levels of aging, immunosenescence, and immune dysfunction biomarkers during treatment was carried out. An improvement with an increase in the level of Natural Killer Cells was reported. After 2 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of Natural Killer Cells improved from 47 to 120.

[00238] After the 2 months of treatment, during which time the patient was administered G-CSF as in Example 18 SPECT Scans (Figure 24) showed improvement in the patients Neurodegenerative Disease.

[00239] After the 2 months of treatment, during which time the patient was administered G-CSF as in Example 18 the patient's insulin requirement for Type 1 Diabetes Mellitus decreased by 50%.

[00240] While not wishing to be bound to any particular theory of operation, increase in Natural Killer Cells and results of SPECT Scans were indicative of improvement in clinical symptoms of immune dysfunction and immunosenescence associated with Chronic disease, Chronic Fatigue Syndrome, Neurodegenerative Disease and Type 1 Diabetes Mellitus.

[00241] Regarding Patient 8, reference is made to Figure 25.

[00242] Diagnosis: Chronic disease, Chronic Viral Infection

[00243] In accordance with the basic protocol outlined in Example 18, a patient BF was treated with stem cell mobilization factor, G-CSF to improve the levels of aging biomarkers in the patient. Evaluation of the improvement in the levels of aging, immunosenescence, and immune dysfunction biomarkers during treatment was carried out. An improvement with an increase in the level of Natural Killer Cells Activity was reported. After 2 months of treatment, during which time the patient was administered G-CSF as in Example 18, levels of Natural Killer Cells Activity improved from 3.02 to 35.44.

[00244] While not wishing to be bound to any particular theory of operation, increase in Natural Killer Cells resulted in improvement in clinical symptoms of immune dysfunction and immunosenescence associated with Chronic disease and Chronic Viral Infection.

[00245] Regarding Patient 9, reference is made to Figures 26-29.

[00246] Diagnosis: Chronic disease, Cancer (Colon Cancer), Frailty.

[00247] In accordance with the basic protocol outlined in Example 5, a patient LK was treated with stem cell activation with stem cell mobilization with granulocyte colony stimulating factor (G-CSF) and in combination with infusing stem-cell rich plasma from ABO-matched healthy allogeneic donors to improve the levels of anti-aging biomarkers in the recipients.

[00248] Evaluation of the improvement in the levels of aging, immunosenescence, immune dysfunction, and early lymphoid lineage differentiation biomarkers during treatment was carried out. An improvement with an increase in the level of the Naive CD4 levels and a decrease in the Memory CD4 T cells and improvement in Natural Killer Cell Activity was reported. After 12 months, during which time the patient was administered G-CSF as in Example 5, levels of Naive CD4 cells improved from 30.00 to 51.79

[00249] Evaluation of the improvement in immune dysfunction resulted in a decrease in the levels of Memory CD4 T Cells from 60.24 to 38.53.

[00250] Also an improvement in Natural Killer Cell Activity was reported After 12 months, during which time the patient was treated as in Example 5, wherein levels of Natural Killer Cell Activity improved from 7.21 to 19.26.

[00251] Also an improvement in B Cell was reported. After 12 months of treatment, during which time the patient was treated as in Example 5, wherein levels of B Cells improved from 36 to 46.

[00252] While not wishing to be bound to any particular theory of operation, reduction in Central Memory T cells and increase in Naive T cells, B Cells and Natural Killer Cell Activity resulted in improvement in clinical symptoms of immune dysfunction , immunosenescence and impairment of early lymphoid differentiation associated with Chronic disease, Cancer (Colon Cancer),Frailty.

[00253] All patents and publications mentioned in this specification are indicative of the levels of those skilled in the art to which the invention pertains. All patents and publications are herein incorporated by reference to the same extent as if each individual publication was specifically and individually indicated to be incorporated by reference.

[00254] It is to be understood that while a certain form of the invention is illustrated, it is not to be limited to the specific form or arrangement herein described and shown. It will be apparent to those skilled in the art that various changes may be made without departing from the scope of the invention and the invention is not to be considered limited to what is shown and described in the specification and any drawings/figures included herein.

[00255] One skilled in the art will readily appreciate that the present invention is well adapted to carry out the objectives and obtain the ends and advantages mentioned, as well as those inherent therein. The embodiments, methods, procedures and techniques described herein are presently representative of the preferred embodiments, are intended to be exemplary and are not intended as limitations on the scope. Changes therein and other uses will occur to those skilled in the art which are encompassed within the spirit of the invention and are defined by the scope of the appended claims. Although the invention has been described in connection with specific preferred embodiments, it should be understood that the invention as claimed should not be unduly limited to such specific embodiments. Indeed, various modifications of the described modes for carrying out the invention which are obvious to those skilled in the art are intended to be within the scope of the following claims.

CLAIMS

What is claimed is:

1. G-CSF for use in treating or preventing a disease associated with aging in a patient, wherein said G-CSF is for administration with an electromagnetic signal and/or a stem cell containing composition.
2. A stem cell containing composition for use in treating or preventing a disease associated with aging in a patient, wherein said stem cell containing composition is for administration with an electromagnetic signal and/or G-CSF.
3. An electromagnetic signal for use in treating or preventing a disease associated with aging in a patient, wherein said electromagnetic signal is for administration with G-CSF and/or a stem cell containing composition.
4. The use of G-CSF in the manufacture of a medicament for treating or preventing a disease associated with aging in a patient, wherein said medicament is for administration with an electromagnetic signal and/or a stem cell containing composition.
5. The use of a stem cell containing composition in the manufacture of a medicament for treating or preventing a disease associated with aging in a patient, wherein said medicament is for administration with an electromagnetic signal and/or G-CSF.
6. The G-CSF for use according to claim 1, the stem cell containing composition for use according to claim 2, the electromagnetic signal for use according to claim 3, the use according to claim 4 or claim 5, wherein said administration is simultaneous, separate, or sequential.
7. A method of treating or preventing a disease associated with aging in a patient, comprising co-administration of G-CSF and an electromagnetic signal to a patient, optionally wherein said method further comprises co-administration of a stem cell containing composition.

8. A method of treating or preventing a disease associated with aging in a patient, comprising co-administration of a stem cell containing composition and an electromagnetic signal to a patient, optionally wherein said method further comprises co-administration of G-CSF.
9. The method according to claim 7 or claim 8, wherein said co-administration is simultaneous, separate, or sequential.
10. A medicament for treating or preventing a disease associates with aging in a patient, comprising G-CSF and a stem cell containing composition.
11. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to any one of the preceding claims, wherein said disease is selected from the group consisting of immunosenescence, immune dysfunction, inflammation, and impairment of early lymphoid lineage differentiation.
12. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to any one of the preceding claims, wherein said disease is selected from the group consisting of anemia, chronic diseases, autoimmune disorders, cancer, cardiovascular diseases, infection, metabolic diseases, neurodegenerative diseases, protein energy malnutrition and frailty.
13. A method of mobilising stem cells and/or stem cell activation in a patient, said method comprising administering to a patient an electromagnetic signal and at least one of:
 - (i) G-CSF; and
 - (ii) a stem cell containing composition.
14. A method of mobilising stem cells and/or stem cell activation in a patient, said method comprising administering to a patient G-CSF and a stem cell containing composition, optionally wherein said method further comprises administering an electromagnetic signal.

15. A method of enhancing the concentration of stem cells at a disease site in a patient, said method comprising administering to a patient an electromagnetic signal and at least one of:

- (i) G-CSF; and
- (ii) a stem cell containing composition.

16. A method of enhancing the concentration of stem cells at a disease site in a patient, said method comprising administering to a patient G-CSF and a stem cell containing composition, optionally wherein said method further comprises administering an electromagnetic signal.

17. A method of reversing, preventing or treating replicative senescence in a patient, said method comprising administering to a patient an electromagnetic signal and at least one of:

- (i) G-CSF; and
- (ii) a stem cell containing composition.

18. A method of reversing, preventing or treating replicative senescence in a patient, said method comprising administering to a patient G-CSF and a stem cell containing composition, optionally wherein said method further comprises administering an electromagnetic signal.

19. The method according to any one of claims 13 to 18, wherein said administering is simultaneous, separate, or sequential.

20. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to any one of the preceding claims, wherein said stem cell containing composition comprises allogenic and/or autologous stem cells.

21. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to claim 20, wherein said composition comprises plasma containing autologous stem cells.

22. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to claim 20, wherein said composition comprises plasma containing allogeneic stem cells.
23. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to claim 22, wherein said allogeneic stem cells are derived from cord blood.
24. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to any one of the preceding claims, wherein said electromagnetic signal is an oscillating magnetic signal.
25. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to any one of the preceding claims, wherein said electromagnetic signal corresponds to a natural magnetic field waveform.
26. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to claim 25, wherein said natural magnetic field waveform is emitted from a patient with a disease status, optionally wherein said disease status and the disease to be treated are the same.
27. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to any one of the preceding claims, wherein said electromagnetic signal is the difference between the natural magnetic field waveform of a healthy individual and the natural magnetic field waveform of a patient with a disease status, optionally wherein said disease status and the disease to be treated are the same.
28. The G-CSF for use, the stem cell containing composition for use, the electromagnetic signal for use, the use, the method or the medicament according to any one of the preceding

claims, wherein said signal is administered or is for administration to a patient by passing a current through a wire loop assembly, optionally wherein said wire loop assembly is selected from the group consisting of a solenoid, a toroid, a planar or flat coil, and a Helmholtz coil.





FIGURE 1

**FIGURE 2**




FIGURE 3





FIGURE 4



Natural complex biologic waveforms of a frog nerve

FIGURE 5





FIGURE 6



Magnetic field generation by current flow

FIGURE 7




FIGURE 8




FIGURE 9




FIGURE 10




FIGURE 11





FIGURE 12

**FIGURE 13**

Patient 2




FIGURE 14

Patient 2




FIGURE 15

Patient 3




FIGURE 16

Patient 4




FIGURE 17

Patient 4




FIGURE 18

Patient 5




FIGURE 19

Patient 5




FIGURE 20

Patient 5




FIGURE 21

Patient 6




FIGURE 22

Patient 7




FIGURE 23

Patient 7




FIGURE 24

Patient 8




FIGURE 25

Patient 9




FIGURE 26

Patient 9




FIGURE 27

PATIENT 9





FIGURE 28

Patient 9

**FIGURE 29**