



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁶ : A23L 1/28, 1/227, 1/229	A1	(11) International Publication Number: WO 99/30579 (43) International Publication Date: 24 June 1999 (24.06.99)
(21) International Application Number: PCT/GB98/03716 (22) International Filing Date: 11 December 1998 (11.12.98) (30) Priority Data: 9726452.7 16 December 1997 (16.12.97) GB (71) Applicant (for all designated States except US): ZENECA LIMITED [GB/GB]; 15 Stanhope Gate, London W1Y 6LN (GB). (72) Inventors; and (75) Inventors/Applicants (for US only): RODGER, Graham, Wood [GB/GB]; Marlow Foods, Station Road, Stokesley TS9 7AB (GB). CORDELL, Geoffrey, Bryan [GB/GB]; Marlow Foods, Belasis Site, Billingham TS23 1YN (GB). MOT-TRAM, Donald, Stewart [GB/GB]; The University of Reading, Dept. of Food Science & Technology, Whiteknights, Reading RG6 6AP (GB). (74) Agents: LOCKE, Timothy, John; Zeneca Specialties, Intellectual Property Group, Belasis Avenue, P.O. Box 2, Billingham TS23 1YN (GB) et al.		(81) Designated States: AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG). Published <i>With international search report.</i> <i>Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.</i>
(54) Title: FLAVOURING MATERIALS FROM YEAST EXTRACTS (57) Abstract Flavouring materials may be produced from filamentous fungi by contacting them with water at a temperature sufficient to reduce their nucleic acid content and concentrating or separating solids from the resulting aqueous solution. The materials may be further subjected to a chemical reaction for example with a sulphur containing amino acid.		

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece			TR	Turkey
BG	Bulgaria	HU	Hungary	ML	Mali	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MN	Mongolia	UA	Ukraine
BR	Brazil	IL	Israel	MR	Mauritania	UG	Uganda
BY	Belarus	IS	Iceland	MW	Malawi	US	United States of America
CA	Canada	IT	Italy	MX	Mexico	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NE	Niger	VN	Viet Nam
CG	Congo	KE	Kenya	NL	Netherlands	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NO	Norway	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	NZ	New Zealand		
CM	Cameroon			PL	Poland		
CN	China	KR	Republic of Korea	PT	Portugal		
CU	Cuba	KZ	Kazakstan	RO	Romania		
CZ	Czech Republic	LC	Saint Lucia	RU	Russian Federation		
DE	Germany	LI	Liechtenstein	SD	Sudan		
DK	Denmark	LK	Sri Lanka	SE	Sweden		
EE	Estonia	LR	Liberia	SG	Singapore		

FLAVOURING MATERIALS FROM YEAST EXTRACTS

THIS INVENTION relates to flavouring materials.

It is known to use hydrolysed yeast extracts as flavouring materials. Yeast is non-toxic to humans and is normally cultured at high density (high dry cell weight per litre).

5 The nucleic acid content of filamentous fungi may be reduced by contacting them with water at high temperatures and separating them from the water, and such a process is described for *Fusarium* in PCT patent Application WO95/23843. We have discovered that the water from which the fungus is separated contains materials which can be used as or converted to flavouring materials for foods, especially if the fungus is *Fusarium*, for example
10 *Fusarium* IMI 145,425.

The current invention comprises a method whereby the soluble components lost from filamentous fungal cells as a result of this heat treatment can be isolated and used as, or converted into, flavouring substances for foods.

15 This invention comprises a method of processing a filamentous fungus to improve its suitability as food which comprises subjecting it in the presence of water to a temperature sufficient to reduce its nucleic acid content substantially characterised by using materials removed from the fungus in the said method directly or after chemical reaction to flavour food.

20 The invention also comprises a flavouring material for food which is an aqueous solution which comprises nucleic acids removed from a filamentous fungus by contacting it with water at an elevated temperature in which the concentration of dissolved solids is sufficiently high to render the material stable to storage at a temperature of 20°C for a period of one month or is a solid comprising such nucleic acids or is a flavouring material comprising a reaction product of such nucleic acids with a sulphur containing amino acid,
25 hydrogen sulphide or ammonium sulphide.

The flavouring materials when in the form of an aqueous solution preferably comprises at least 30% by weight and more preferably 45 to 60% by weight of solids.

Whilst taste is an important factor in food flavours, the odours of flavouring materials are also important.

30 The soluble components are preferably concentrated from the aqueous solution arising from the nucleic acid reduction step by removing water for example by evaporation, distillation (preferably at reduced pressure) reverse osmosis, freeze drying or freezing out the water as ice leaving an aqueous concentrate. It may suitably be removed by evaporation at reduced pressure for example at a temperature of 40 to 70°C

35 The dissolved solids may be separated as such or left as a concentrated solution where the *A_w* (water activity) is reduced sufficiently to ensure biostasis at a range of storage temperatures.

If the nucleic acid content of the filamentous fungus is reduced by raising the temperature of its growth medium the water recovered will contain salts and other nutrients,

for example glucose and/or complex nitrogen nutrients in addition to the nucleic acids and other materials derived from the fungus. If the flavour imparted by such materials is required they may be left in the materials, but if not they may either be removed, for example by osmosis or ultrafiltration, or the fungus may be washed before its nucleic acid content is reduced thereby avoiding their presence. In WO95/23843 the removal of nucleic acid from a filamentous fungus in its growing state is described; such a process is an improvement over the treatment of fungus in its resting state, for example in pure water. We have found however that the organism takes a short time to adjust from its growing to its resting state and that providing the nucleic acids are removed soon after it is separated from its growth medium the nucleic acids may be satisfactorily removed according to the procedure of WO95/23843.

We have found that after partial or complete removal of water as aforesaid the concentrate can be used as an alternative to hydrolysed vegetable proteins, yeast autolysates or yeast extract as an additive for food. The materials removed from the fungus of value in the production of savoury flavouring preparations and process flavourings. Because of the savoury nature of the flavour it may be used directly in the flavouring of Snacks, Biscuits, Stocks, Soups, Stews, Sauces and Gravies at inclusion levels of preferably between 0.1 and 15 for example 1 to 10 dry weight %.

We have also found that on heating it produces an attractive roast-type aroma. If desired it may be partially hydrolysed before heating, for example by hydrolysis with acetic acid, to produce modified roast flavours.

It may also be reacted, optionally after at least partial hydrolysis, with sulphur containing amino acids, preferably cysteine or optionally with H_2S and/or $(\text{NH}_4)_2\text{S}$ to produce savoury flavours.

The savoury nature of the material may be altered by chemical reaction to provide a different flavour profile in that the meaty/roasted flavour notes are increased. Such "reaction flavourings" may be used in flavouring Meat (beef, chicken, lamb, pork, etc), meat alternatives (e.g. based on soy, wheat, pea protein, myco-protein), prepared meals, snacks and drinks at inclusion levels of preferably between 0.1 and 10 for example 1 to 8 dry weight %.

The flavourings may be produced by reacting materials removed from the fungus as aforesaid with cysteine. This may be carried out in the presence of water if desired; for example a 1.5 to 75 and preferably a 5 to 50 weight % solution of such materials may be reacted with cysteine in quantities of up to 10%, for example 1 to 5% of cysteine by weight based on such materials. The reaction may be carried out at a temperature of for example 110 to 140°C at a pH of 5.5 to 9. Reaction is suitably continued for 0.5 to 7.5 hours.

It is believed that hydrolysis increases the free ribose content of the concentrate and this may be appropriate if certain flavours are desired. It is desirable to avoid treatment with

hydrochloric acid for regulatory reasons (possible production of chloro-propanol derivatives), but hydrolysis with for example acetic acid may be desirable.

In the following descriptions the term "Centrate" is used for the extracellular liquid recovered by the heat shock treatment of a suspension of *Fusarium* at about 70°C in the presence of its growth medium after separation of the cellular material. The term FDC means "freeze dried centrate".

EXPERIMENTAL PROCEDURES

Material Preparation

The liquid centrate was freeze dried in order to:

- reduce the water content and therefore inhibit microbial growth;
- carry out studies at a high concentration of centrate;
- facilitate the handling of the product.

All further analyses described in this report deal with the freeze dried centrate abbreviated FDC.

Methods - Compositional Analyses

Moisture

The moisture content was determined by measuring the weight decrease of the FDC, until constant weight, while placed in an oven at 100°C.

Ash

The ash content was determined by placing the FDC in an oven at 600°C until constant weight was obtained.

Organic Nitrogen

Kjeldahl nitrogen determination was carried out; sucrose was used as blank and glycine as standard. The results are shown in Table 1.

Table 1 - Moisture, Ash and Organic Nitrogen Content of FDC

	Replicates			Mean
	1	2	3	
Moisture (%)	13	13	13	13
Ash (%)	18	18	17	18
Organic Nitrogen (%)	6	6	6	6

Amino Acids

The amino acids determination was done with a 6300 Beckman auto-analyser. The free amino acids present in the FDC were analysed using a 0.06% solution of the FDC. It was possible to measure the total amino acid content by a prior hydrolysis (HCl 6N, 24h, oven 110°C) of the FDC. However, the acid is known to hydrolyse tryptophan and the sulphur amino acids. The hydrolysis of the sulphur amino acids can be avoided by a prior oxidation of cysteine into cysteic acid, and methionine into methionine sulfone. This was

carried out by treating the FDC with a solution of formic acid/hydrogen peroxide/methanol (48.5/1/0.5) during 4h at 0°C in the dark. The results are shown in Table 1a.

Table 1a - Amino Acid Content of FDC

mean results in g/100g FDC

Amino Acids	FDC (free AA)	Hydrolysed FDC (total AA)
CYSTEIC ACID	0.09	0.32
ASP	0.14	0.50
THR	0.02	0.24
SER	0.12	0.16
GLU	2.02	1.70
CYSTEINE	0.00	0.21
PRO	0.00	0.20
GLY	0.05	0.24
ALA	2.05	1.73
VAL	0.14	0.37
CYSTINE	0.00	0.01
METH	0.00	0.13
ILE	0.00	0.21
LEU	0.00	0.30
TYR	0.00	0.06
PHE	0.00	0.17
TRYPTOPHAN	0.00	0.00
NH ₃	1.49	4.56
LYS	0.15	0.26
HIS	0.00	0.08
ARG	0.71	0.65
TOTAL	6.98	11.80

5 **Total Carbohydrates**

The carbohydrate content of the FDC was assessed by the phenol-sulphuric acid assay method (Carbohydrate analysis: a practical approach, ed. Chaplin, Kennedy, IRL Press). Solutions of FDC and glucose (standard for calibration) were mixed with a solution of phenol in water (5% w/v). Concentrated sulphuric acid (1 ml) was added rapidly and directly to the solution surface without allowing it to touch the sides of the tube. The solutions were left undisturbed for 10 min before shaking vigorously. The absorbencies were read at 490 nm after a further 30 min.

Sugars

The sugar analysis was performed using a Dionex System of High Pressure Liquid Chromatography (HPLC), in which an eluent of HPLC grade water (1 ml/min) was used with an anion-exchange column (column Dionex PA-1) and a pulsed amperometric detector.

Pure compounds were used as standards for retention time determination and quantitation.

Free sugars were analysed using a 0.15% solution of FDC after filtering the solution through a 0.45 µm Minisart 25 membrane. Total sugars were also evaluated after a preliminary acid hydrolysis of the FDC (solution 0.15% in HCl 1N, 2h, oven 110°C) and filtration through first an Ag filter (precipitate of AgCl) and second a 0.45 µm Minisart 25 filter. The results are shown in Table 2.

Table 2 - HPLC Analysis of Sugar Content of FDC

g sugar per 100g of FDC

	FDC mean	HYDROLYSED FDC mean
arabinose	0.02	0.04
galactose	0.10	0.10
glucose	8.62	23.01
sucrose	0.09	0.23
xylose	0.00	0.00
mannose	0.00	0.10
fructose	0.01	0.00
ribose	0.00	0.06
maltose	2.08	
TOTAL SUGARS	12.66	25.10

Nucleic acid Derivatives

A Perkin Elmer Binary HPLC pump 250 equipped with a Spectroflow 757 ABI Analytical Kratos Division was used. Standards and samples were filtered through Acrodisc 0.45 µm Gelman Sciences membranes filters and injected by means of an injector valve equipped with a 20 µl injection loop into a reverse-phase µBondapack C18 (3.9 x 300 mm) Waters analytical column, protected by a µBondapack C18 guard column. A wavelength of 254 nm was used. A gradient programme with two mobile phases was used: mobile solvent A was a 60/40 methanol/water mixture and mobile solvent B was 0.02M KH₂PO₄ (pH 5.5) prepared from potassium dihydrogen orthophosphate in distilled water and pH adjusted with 1M KOH. All mobile solvents were filtered (Nylaflo 0.2 µm Gelman Sciences membranes filters) and degassed with Helium before use. The total run time was 51 min and the flow rate 1 ml/min which consisted of 100% solvent B during 5 min, followed by a gradient from 0% to 36% solvent A in 36 min and 36% solvent A for 5 min. Then a reverse

gradient of 36% to 0% A was set for 5 min and the HPLC was ready for further injection after 15 min equilibrium.

Identification of the compounds was made by comparison with the retention time obtained from standards analysed in the same HPLC conditions. Standards were analysed separately to know their individual retention time and then all together to check any elution over lap that may occur in the sample case. These standards are presented in Table 3.

Hydrolysis of FDC

Since free ribose is highly reactive in the Maillard reaction, the effect of gentle hydrolysis conditions on FDC and the subsequent effects on flavour generation were investigated. Acid hydrolysis was carried out with sodium acetate 0.01M, pH4 adjusted with acetic acid. Standards (inosine, adenosine 5' mono phosphate -AMP5', guanosine and guanosine 5' mono phosphate -GMP5') were prepared at 4000 μ M in duplicate and an aliquot of each solution was taken and run under the same HPLC conditions as were adapted for the analysis of nucleic acids derivatives above. The solutions were then subjected to hydrolysis for 7.5 h in an oven 100°C (GC oven Carlo Erba). The reaction was stopped by placing the tubes in an ice bath and kept in freezer until analysis.

FDC at 2% w/v was subjected to similar hydrolysis conditions.

Flavour mixture preparation

As indicated previously, FDC is considered to have the potential of being either a flavouring in its own right, or a precursor in the generation of reaction - product flavours. Therefore a range of reaction mixtures were prepared and are presented in Table 4.

Table 4 - Flavoured mixtures analysed by sniffing panel
 Selected mixtures for GS-MS analysis are underlined
 C = Cysteine

Sample Name Sample Composition	<u>Aqueous Centrate</u> 1.7%(w/v, solids of centrate/water) 0.5h, 140°C	<u>Heated Aqueous Centrate</u> 1.7%(w/v, solids of centrate/water) 0.5h, 140°C	
Sample Name Sample Composition	<u>Heated Buffered (pH 5.5) Centrate</u> 1.7% (w/v, solids of centrate/sodium acetate 0.01M) 0.5h, 140°C	<u>Heated Hydrolysed Buffered (pH 5.5) Centrate</u> 1.7% (w/v, solids of centrate/sodium acetate 0.01M) 7.5h, 110°C 0.5h, 140°C	<u>Heated Hydrolysed Buffered (pH 5.5) Centrate + C</u> 1.7% (w/v, solids of centrate/sodium acetate 0.01M) 1g cysteine/17g solids of centrate 7.5h, 110°C 0.5h, 140°C
Sample Name Sample Composition	<u>Heated Aqueous Centrate 12% pH 5.5</u> 12% (w/v, solids of centrate/water) 0.5h, 140°C	<u>Heated Aqueous Centrate 20% pH 5.5</u> 20% (w/v, solids of centrate/water) 0.5h, 140°C	<u>Heated Aqueous Centrate 30% pH 5.5</u> 30% (w/v, solids of centrate/water) 0.5h, 140°C
Sample Name Sample Composition	<u>Heated Aqueous Centrate 50% pH 5.5</u> 50% (w/v, solids of centrate/water) 0.5h, 140°C	<u>Heated Aqueous Centrate 75% pH 5.5</u> 75% (w/v, solids of centrate/water) 0.5h, 140°C	<u>Heated Aqueous Centrate 87% pH 5.5</u> 87% (w/v, solids of centrate/water) 0.5h, 140°C
Sample Name Sample Composition	<u>Heated Aqueous Centrate 20% pH 7.5</u> 20% (w/v, solids of centrate/water) 0.5h, 140°C	<u>Heated Aqueous Centrate 30% pH 7.5</u> 30% (w/v, solids of centrate/water) 0.5h, 140°C	
Sample Name Sample Composition	<u>Heated Aqueous Centrate 20% pH 9</u> 20% (w/v, solids of centrate/water) 0.5h, 140°C	<u>Heated Aqueous Centrate 30% pH 9</u> 30% (w/v, solids of centrate/water) 0.5h, 140°C	
Sample Name Sample Composition	<u>Heated Aqueous Centrate pH 5.5 + C 1/20</u> 20% (w/v, solids of centrate/water)	<u>Heated Aqueous Centrate pH 5.5 + C 1/10</u> 20% (w/v, solids of centrate/water)	<u>Heated Aqueous Centrate pH 5.5 + C 1/5</u> 20% (w/v, solids of centrate/water)

	1g cysteine/20g solids of centrate 0.5h, 140°C	1g cysteine/10g solids of centrate 0.5h, 140°C	1g cysteine/5g solids of centrate 0.5h, 140°C
Sample Name Sample Composition	<u>Heated Aqueous Centrate pH 9 + C</u> <u>1/20</u> 20% (w/v, solids of centrate/water) 1g cysteine/20g solids of centrate 0.5h, 140°C	<u>Heated Aqueous Centrate pH 9 + C</u> <u>1/10</u> 20% (w/v, solids of centrate/water) 1g cysteine/10g solids of centrate 0.5h, 140°C	<u>Heated Aqueous Centrate pH 9 + C</u> <u>1/5</u> 20% (w/v, solids of centrate/water) 1g cysteine/5g solids of centrate 0.5h, 140°C
Sample Name Sample Composition	<u>175°C Heated Aqueous Centrate</u> <u>pH 5.5 + C 1/20</u> 20% (w/v, solids of centrate/water) 1g cysteine/20g solids of centrate 5 min, 175°C	<u>100°C 1h Heated Aqueous</u> <u>Centrate pH 5.5 + C 1/20</u> 20% (w/v, solids of centrate/water) 1g cysteine/20g solids of centrate 1h, 100°C	<u>100°C 1.5h Heated Aqueous</u> <u>Centrate pH 5.5 + C 1/20</u> 20% (w/v, solids of centrate/water) 1g cysteine/20g solids of centrate 1.5h, 100°C
Sample Name Sample Composition	<u>175°C Heated Aqueous Centrate</u> <u>pH 9 + C 1/5</u> 20% (w/v, solids of centrate/water) 1g cysteine/5g solids of centrate 5 min, 175°C	<u>100°C 1h Heated Aqueous</u> <u>Centrate pH 9 + C 1/5</u> 20% (w/v, solids of centrate/water) 1g cysteine/5g solids of centrate 1h, 100°C	<u>100°C 1.5h Heated Aqueous</u> <u>Centrate pH 9 + C 1/5</u> 20% (w/v, solids of centrate/water) 1g cysteine/5g solids of centrate 1.5h, 100°C

Reaction mixtures (2 ml) were prepared by mixing appropriate quantities of stock solutions in glass tubes and then transferring to 20 ml Kimble ampules that were sealed in hot flame. The ampules were then placed in a metal cover and heated in a Carlo Erba 4200 gas chromatograph oven.

5 The reaction mixtures were stored in the freezer at -20°C before analysis. The ampules were broken for analysis after bringing the reaction mixtures to room temperature.

Sensory evaluation of Aroma volatiles

An informal panel of 6 persons (3 females, 3 males) experienced in flavour evaluation was recruited.

10 For sensory evaluation, 1 ml aliquots of the samples under investigation were transferred into brown screw-cap bottles and diluted 10 times (except for the concentration study where no dilution was applied). The coded samples were presented to one panellist at one time at room temperature and the panellists were asked to describe the aromas using their own terms.

Instrumental evaluation of Aroma volatiles determination

15 A dynamic headspace collection procedure was used. Each sample (1.7 ml of reaction mixture so that it was equivalent to 0.4 g of FDC) was placed in a 250 ml conical flask fitted with a Drechsel head. Distilled water was added to a final volume of 10 ml and the mixture shaken gently. Oxygen-free nitrogen was passed over the sample for 1 h at a rate of 40 ml/min. The volatiles were swept onto a preconditioned glass-lined stainless-steel trap (105 mm x 3 mm i.d.) packed with 85 mg Tenax GC (CHIS system, SGE Limited). Throughout the collection, the sample was maintained at 37°C using a water bath. The internal standard was 1,2-dichlorobenzene in ether (130 µl/ml) and 1 µl was injected onto the trap at the end of the collection time, the trap was then flushed with nitrogen for 10 min.

25 A Hewlett-Packard (HP) 5890/5972 gas chromatograph-mass spectrometer (GC-MS), fitted with a 50 m x 0.32 mm i.d. fused-silica capillary column coated with BPX-5 (SGE Limited) at 0.5 µm film thickness, was used to analyse the collected volatiles. These were thermally desorbed at 250°C in the CHIS injection port (SGE Limited) and cryofocused directly onto the front of the GC column, while the oven was held at 0°C for 5 min. The oven temperature was then raised to 40°C over 1 min and held for 5 min before raising the temperature to 250°C at a rate of 4°C/min and holding for a further 10 min. The helium carrier gas flow rate was 1.5 ml/min. Mass spectra were recorded in the electron impact mode at an ionisation voltage of 70 eV and source temperature of 200°C. A scan range of 29-400 m/z and a scan time of 0.69 s were used. The data were controlled and stored by the HP G1034C Chemstation data system.

35 Volatiles were identified by comparison of their mass spectra with the spectra from authentic compounds in the Reading Laboratory or in the NIST/EPA/MSDC Mass Spectral Database or other published spectra. The linear retention index (LRI) was calculated for each component using the retention times of a homologous series of C₆-C₂₂ n-alkanes.

Nucleic Acid Derivatives

The nucleic acid composition of the centrate was determined from 3 replicates and is presented in Table 5. It explains most of the HPLC eluted peaks. A number of compounds co-eluted, but it was not possible to investigate alternative analysis conditions and, therefore, for co-eluting compounds it was not possible to determine which of the compounds contributed to the peak obtained in the centrate analysis.

As expected, there were few deoxyribonucleic acid derived compounds compared with the ribonucleic acid derived compounds, which are more abundant in nature. The major nucleic acid components are cytosine 5' monophosphate (26% of the total nucleic acid content), uridine 3' monophosphate and/or guanosine 5' monophosphate (18%), adenosine 5' monophosphate and/or deoxyribo guanosine 5' monophosphate (16%). All of them are potential sources of ribose and ribose phosphate which are good reactive precursors the Maillard reaction. Excluding the bases, the potential source of ribose or ribose phosphate represents 96% of the nucleic acid content of the centrate, which is equivalent to 202 ppm of the content of the centrate.

Table 5 - Centrate Nucleic Acids HPLC Analysis Results

Nucleic Acid Type	CMP'5	UMP5'	Uracil'	UMP3' + GMP5'	IMP5'	Cytidine	Hypox+ guanine + uridine + GMP3'	AMP5' + I + DAMP3' + DAMP5'	Purine	AMP5' + I + DAMP3' + DAMP5'	DG	DI	Adeno-sine	Total
nucleic acid (ppm in centrate)	54	23	1	38	5	5	7	34	2	6	2	15	19	212
nucleic acid (% of total)	26	11	0	18	3	2	3	16	1	3	1	7	9	100
Standard deviation (ppm in centrate)	3	1	0	1	0	1	0	2	0	0	0	1	2	

hypox = hypoxanthine

I = inosine

C = cytidine

U = Uridine

G = Guanosine

A = Adenosine

MO = mono phosphate

D = deoxyribose

Effects of Acid Hydrolysis on Nucleic Acid Derivatives

Table 6 presents the results of hydrolysis of solutions of inosine, guanosine and their respective 5' phosphate ribonucleotides. The method is based on that used by Matoba et al (J. Food Science, vol. 53, n.4, 1988, p1156). The last column gives an indication of quantity recovery and it can be seen that the results of the hydrolysis on the

The most interesting model systems are the ribonucleotides since they are major components in the centrate. They were hydrolysed by half or less, producing their respective nucleosides which were further hydrolysed into their bases. Although it is possible to hydrolyse the ribonucleotides into their bases and consequently produce ribose and/or ribose phosphate, relatively low yields of bases were obtained and an optimisation of this process should be carried out.

Table 6 - Acid hydrolysis results of model systems

	Quantity (%)	After hydrolysis			
	initial	hypoxanthine	inosine	IMP5'	total
inosine	100	13	85		98
IMP5'	100	6	15	71	92
		guanine	guanosine	GMP5'	
guanosine	100	2	42		44
GMP5'	100	2	35	51	88

Conclusion

The nucleic acid composition of the FDC has been characterised. It comprises mainly ribonucleotides with relatively small amounts of deoxyribonucleotides. Hydrolysis of nucleotides releases free ribose or ribose phosphate only occurs to a relatively small extent in acetate buffer at pH 4.

Results - Sensory Evaluations of Aroma Volatiles

It was decided to present the individual results of each panellist and not to group them under specific common descriptors because of the too large diversity of the terms described.

Tables 7 and 8 present the effect of heating and the impact of the hydrolysis, with or without the addition of cysteine.

Table 7 - Aroma panel results on centrate - Study of the effect of cooking

Sample name	Aqueous Centrate	Heated Aqueous Centrate
Sample Composition	1.7%(w/v, solids of centrate/water	1/7%(w/v, solids of centrate/water) 0.5, 140°C
Panellist 1	scrumpy, glucose, syrup, molasses	molasses, caramel
Panellist 2	caramel, sweet	burnt
Panellist 3	whey, old yoghurt, sharp, creamy	raw celery, braised celery, weird smell
Panellist 4	wet cloth/ironing/scorching auto-claving media, slightly acrid	fermenting cereal, ironing/wet cloth, cotton/wool, treacle, golden syrup
Panellist 5	honey, urine	caramel, fatty
Panellist 6	caramel, slightly fruity	slightly fruity caramel, burnt, sharp

Table 8 - Aroma panel results on heated buffered centrate - Study of the effect of hydrolysis and of addition of cysteine

Sample	Heated Buffered (pH 5.5) Centrate	Heated Hydrolysed Buffered (pH 5.5) Centrate	Heated Hydrolysed Buffered (pH 5.5) Centrate + Cysteine
Sample Composition	1.7% (w/v) solids of centrate/sodium acetate 0.01M) 0.5h, 140°C	1.7% (w/v, solids of centrate/sodium acetate 0.01M) 7.5h, 110°C 0.5h, 140°C	1.7% (w/v, solids of centrate/sodium acetate 0.01M) 1g cysteine/17g solids of centrate 7.5h, 110°C 0.5h, 140°C
Panelist 1	burnt, caramel, toffee, acid, acid	resinous, burnt, acrid	burning paper/plastic/hair, putrid, acrid
Panelist 2	burnt, caramel, celery, slightly sweet	stale shall, burnt, acidic	burnt, celery, acid
Panelist 3	braised celery, marmite, jammy, cooked apple	varnish, paint, biscuit	sweet, grassy, herbal, raw onion, vinegar
Panelist 4	rotting vegetable, sulfur, autoclaving malted barley, burnt, nutty, treacle	wet cloth, autoclaving, fermenting cereal, nose catching/sharp	onion rotting, hard boiled egg, malty/barley/fermenting cereal
Panelist 5	treacle, honey, fatty	nicotine	burning tyres
Panelist 6	caramel, toffee, slight fruity, slightly burnt	chemical, slightly burnt rubbery	caramel, toffee, sweet, floral

The study of the effect of concentration of centrate involved the range of concentrations likely to be reached in commercial practice, i.e., within the range 12% and 30% of solids. These were compared with the non diluted freeze dried centrate powder (87% solids). The other preparation conditions were kept constant. viz. pH 5.5 and heating at 140°C for 30 min. The results are presented in Table 9. The odours were very strong and the reproducibility of the results within 2 replicates for each panellist was fairly poor. However, there was a noticeable trend within the sample set from low to high concentration: at 12%, the odours were mainly sweet, vegetable and molasses. These notes became associated with burnt and sharp as well as savoury at 20 and 30% solids. At 50%, the sample had roasted and paint smells that became dominant at 75%. Some extra metallic, burnt rubbery, and sulphur notes were detected with the 87% sample. The 20 and 30% solid samples seemed the most interesting because of their meaty savour smells and therefore they were selected for further analysis. it was also decided to dilute the original flavoured reaction mixtures before sniffing further reaction mixtures.

Table 9 - Aroma panel results on heated concentrate pH 5.5 - Study of the effect of concentration

Sample Name	Heated Aqueous Concentrate 12% pH 5.5	Heated Aqueous Concentrate 20% pH 5.5	Heated Aqueous Concentrate 30% pH 5.5
Sample Composition	12% (w/v, solids of concentrate/water) 0.5h, 140°C	20% (w/v, solids of concentrate/water) 0.5h, 140°C	30% (w/v, solids of concentrate/water) 0.5h, 140°C
Panellist 1	burnt, caramel, slightly sweet	sweet, slightly meaty, nutty	soja sauce, caramel, toffee, burnt, sweet
Panellist 2	sickly molasses	sharp, acid, cloying, burnt paper	molasses, black treacle
Panellist 3	celery, burnt roast coffee, bovril	marmite	marmite, emulsion paint
Panellist 4	honey, caramel, stir fry	honey, caramel, stir fry, vegetable	caramelised vegetables
Panellist 5	celery, rancid, maggi, bovril	sharp, marmite, yeasty, bovril+celery undertones	sharp, bovril, vinegary top note
Panellist 6	malty, brewery, virol (malt extract), v. concentrate digestive biscuit, black treacle, autoclaving	honey, malty	autoclaving, honey, digestive biscuits
	REPLICATE	REPLICATE	REPLICATE
Panellist 1	sweet, caramel, soja sauce	burnt, toffee, acid	soja sauce, caramel, burnt, sweet
Panellist 2	burnt, caramel, burnt skin/hair	sharp, acid, cloying, burnt paper	molasses, resinous
Panellist 3	malt	marmite, celery	marmite, celery, slight roast coffee
Panellist 4	honey, caramel, stir fry	honey, earthy, uncooked potato	soy sauce, wood smoke
Panellist 5	celery, meaty	burnt celery	bovril, strange top note
Panellist 6	black treacle, virol, autoclaving	concentrated honey, biscuity digestives, autoclaving	biscuits, honey, green/fruity, green bananas

Sample Name	Heated Aqueous Concentrate 50% pH 5.5	Heated Aqueous Concentrate 7.5% pH 5.5	Heated Aqueous Concentrate 87% pH 5.5
Sample Composition	50% (w/v, solids of concentrate/water) 0.5h, 140°C	75% (w/v, solids of concentrate/water) 0.5h, 140°C	87% (w/v, solids of concentrate/water) 0.5h, 140°C
Panellist 1	burnt, caramel, bitter chocolate	burnt, very strong, chocolate	burnt sugar
Panellist 2	molasses, black treacle, solvent	burnt, solvent, sickly	savoury, burnt, burnt flesh/skin
Panellist 3	paint, marmite, slight celery	burnt	slightly sweet, burnt, marmite
Panellist 4	raw vegetables, caramel	caramel, roasted slight treacle	soy sauce

Panellist 5	strange top note of sharp green/ fruity note/household or dry smell	very burnt sugar + top note	burnt rubber, sulphur, very burnt sugar
Panellist 6	green banana, honey, slight autoclaving, slight virol	burnt milk, burnt sugar (not caramel - really burnt)	dry, dusty sensation, like something burnt black in oven, sharp, faint marmite, black treacle, slightly stale, metallic, rusty steel
	REPLICATE	REPLICATE	REPLICATE
Panellist 1	burnt, caramel	burnt, slightly sweet, caramel	burnt sugar
Panellist 2	burnt, black treacle, slightly acrid	very strong marmite, meat extract, burnt, caramel	molasses, burnt, savoury
Panellist 3	paint, marmite, cereal	marmite, burnt	burnt, marmite
Panellist 4	vegetable, roasted, woodland	damp wood, syrup, honey	vegetable
Panellist 5	strong strange top note, bovril	strong strange top note, burnt sugar	marmite, yeasty, bovril, very concentrated
Panellist 6	black treacle, acrid/burnt honey, virol/malt	fresh sensation, honey, malt	same as other replicate

5 The pH effect was studied on the centrate at 3 different values: 5.5, 7.5 and 9, with the heating conditions kept at 140°C for 30 min. The results are in Table 10. There was not much difference in the results between the 20% and the 30% solids samples. The results at pH 7.5 were similar to those at pH 5.5 and the smells were mainly autoclave and caramel. The odour became burnt with pH 9. Therefore it was decided to carry on the sniffing experiments by selecting the two extreme pHs and to keep only one concentration (20% solids).

Table 10 - Aroma panel results on heated centrate study of the effect of pH

Sample Name	Heated aqueous centrate 20% pH 5.5	Heated aqueous centrate 30% pH 5.5
Sample composition	20%(w/v, solids of centrate/water) 0.5h 140°C	30% (w/v, solids of centrate/water) 0.5h 140°C
Dilution before sniffing	50 for panellist 1 10 for panellist 2	50 for panellist 1 10 for panellist 2
Panellist 1	1-slight autoclaving 2-diacetyl then going to caramel, butter scotch, slightly nutty	1-slightly autoclave, "catching" in nose. 2-then quite a lot of caramel
panellist 2	caramel, butter-like	caramel, sweet
Sample name	Heated aqueous centrate 20% pH 7.5	Heated aqueous centrate 30% pH 7.5
Sample composition	20%(w/v, solids of centrate/water 0.5h 140°C	30%(w/v, solids of centrate/water) 0.5h 140°C
Dilution before sniffing	50	50
Panellist 1	similar to sample 20% pH 5.5 in the way it changes, bit more autoclaving, caramel but not really butterscotch	1-slightly autoclave 2-caramel
Sample name	Heated aqueous centrate 20% pH 9	Heated aqueous centrate 30% pH9
Sample composition	20% (w/v, solids of centrate/water) 0.5h 140°C	30% (w/v, solids of centrate/water) 0.5h 140°C
Dilution before sniffing	10	10
Panellist 2	burnt, baked, roasted, cereals	slightly burnt, caramel sweet

The effect of addition of cysteine was studied on the centrate in solution at 20% solids at pH 5.5 and 9. There were three concentrations of cysteine tested: ratios 1/20, 1/10 and 1/5 of cysteine (g)/centrate solids (g). The heating conditions were kept the same as previously (104°C for 30 min). The results are presented in Table 11. Within the pH 5.5 sample series, the low concentration of cysteine sample lead to a somewhat pleasant odour of sweet, greasy, meaty sauce, that was progressively replaced by roasted and rubber notes as the cysteine concentration increased. At pH 9, the burnt roasted cereals notes already mentioned in the previous experiment were present again with cysteine at low concentration. When the cysteine content increased, the odour became strong and more nutty and then close to savoury, meaty stock. It was therefore decided to select pH 5.5 with cysteine 1/20 sample and pH9 with cysteine 1/5 sample for the next set of experiments.

Table 11 - Aroma panel results on heated centrate pH 5.5 and pH 9 - Study of the effect of addition of cysteine

Sample Name	Heated Aqueous Centrate pH 5.5 + cysteine 1/20	Heated Aqueous Centrate pH 5.5 + cysteine 1/10	Heated Aqueous Centrate pH 5.5 + cysteine 1/5
Sample Composition	20% (w/v, solids of centrate/water) 1g cysteine/20 g solids of centrate 0.5h, 140°C	20% (w/v, solids of centrate/water) 1g cysteine/10g solids of centrate 0.5h, 140°C	20% (w/v, solids of centrate/water) 1g cysteine/5g solids of centrate 0.5h, 140°C
Dilution before Sniffing	10	10	10
Panelist 1	celery, cooked vegetables, burnt wool, sweaty socks	bacon fat, crispy duck (chinese style), sl.sweaty, roasted	sweet, caramel, roasted
Panelist 2	celery (strong), sl. burnt coffee, treacle toffee (sl.)	puffed wheat, sweet+sour	celery, sweet+sour, puffed wheat (strong)
Panelist 3	strong chicken stock, greasy, acidic, rubber	burning rubber, savoury, marmite, malt extract	weak rubbery, acrid
Panelist 4	strong celery (main odour) + other notes - meaty, pork and apple, sweet and sour sauces	medicinal smell, also watery stock (weak) faint celery/sweet'n sour	(sweet and sour), rubber, magi
Panelist 5	black treacle, buttery/caramel, burnt sugar, burnt fried boiled cabbage, marsala, sweet, nutty	buttery, golden syrup, something burnt? (wet wood or toasted nuts, almond, hazelnut)	golden syrup, caramel, toasted nuts (not burnt) sl. sulphur.

Sample Name	Heated Aqueous Concentrate pH 9 + cysteine 1/20	Heated Aqueous Concentrate pH 9 + cysteine 1/10	Heated Aqueous Concentrate pH9 + cysteine 1/5
Sample Composition	20% (w/v, solids of concentrate/water) 1g cysteine/20g solids of concentrate 0.5h, 140°C	20% (w/v, solids of concentrate/water) 1g cysteine/10g solids of concentrate 0.5h, 140°C	20% (w/v, solids of concentrate/water) 1g cysteine/5g solids of concentrate 0.5h, 140°C
Dilution before Sniffing	10	10	10
Panelist 1	burnt, roasted, nutty	sweet caramel, roasted, nutty, burnt	meaty, pork crackling, nutty
Panelist 2	sweet (strong), toffee, biscuity, marzipan	celery (strong), sl. sweet, biscuity	meaty (strong)
Panelist 3	savoury, burnt skin, stock	burning rubber, strong savoury, malt extract	strong chicken stock, chicken fat, sl. acid
Panelist 4	v. strong weird note, stale biscuit crumbs, Farleys rusks, faint sweet'n sour note	fruity, sweet, burnt sugar, vinegar	sulphury, weak stock, nasty note, burnt meat?
Panelist 5	burnt, scorched/burnt wet wool, golden syrup	same as sample 1/20 but stronger burnt, less golden syrup	meaty (beef gravy?), burnt, burnt sugar, plus caramel

The results of the temperature/duration heating conditions study are presented in Table 12. They involve the two selected samples previously described which were then cooked at a lower temperature and longer time: 100°C for 60 min and 90 min, and at a higher temperature but for a shorter time: 175°C for 5 min. Compared to the original heating conditions, similar results were obtained for the pH 5.5 cysteine 1/20 sample by heating 100°C for 60 min. Longer time of heating at 100°C resulted in more roasted burnt notes and a similar type of odour was obtained after 5 min at 175°C. Regarding the pH 9 cysteine 1/5 sample, the meaty notes obtained by 30 min at 140°C were reached by the treatment 5 min at 175°C. At 100°C, some very strong odours of urine and boiled eggs were described.

**Table 12 - Aroma panel results on heated centrate pH 5.5 and pH 9 + Cysteine
Study of the effect of temperature and duration of the cooking**

Sample Name	175°C Heated Aqueous Centrate pH 5.5 + cysteine 1/20	100°C 1h Heated Aqueous Centrate pH 5.5 + cysteine 1/20	100°C 1.5h Heated Aqueous Centrate pH 5.5 + cysteine 1/20
Sample Composition	20% (w/v, solids of centrate/water) 1g cysteine/20g solids of centrate) 5min, 175°C	20% (w/v, solids of centrate/water) 1g cysteine/20g solids of centrate 1h, 100°C	20% (w/v, solids of centrate/water) 1g cysteine/20g solids of centrate 1.5h, 100°C
Dilution before sniffing	10	10	10
Panelist 1	burnt-strong, caramel-moderate, roasted-moderate	fruity-moderate, nutty-slight	celery-moderate, roasted-light
Panelist 2	puffed wheat-strong, marmite- medium, celery-light	celery-strong	celery-medium, cereal-slight
Panelist 3	rancid-v. strong, acrid-strong, burnt fat/skin-medium, molasses- medium/weak	chicken fat-medium, rubber- medium, burnt rubber-medium	burnt rubber-weak, molasses-weak medium, burnt/cold wood-weak

Sample Name	175°C Heated Aqueous Centrate pH 9 + cysteine 1/5	100°C 1h Heated Aqueous Centrate pH 9 + cysteine 1/5	100°C 1.5h Heated Aqueous Centrate pH 9 + cysteine 1/5
Sample Composition	20% (w/v, solids of centrate/water) 1g cysteine/5g solids of centrate 5min, 175°C	20% (w/v, solids of centrate/water) 1g cysteine/5g solids of centrate 1h, 100°C	20% (w/v, solids of centrate/water) 1g cysteine/5g solids of centrate 1.5h, 100°C
Dilution before sniffing	10	10	10
Panelist 1	crackers-moderate, stale-slight, baked bread-moderate	yeast extract-weak, vegetables- moderate, wet washing-moderate	yeast, bread dough-moderate, vegetables-weak, nutty-background after other odours decrease
Panelist 2	sulphury-slight, meaty-medium, burnt-medium	hard boiled egg-strong	egg-strong
Panelist 3	burnt paper-weak, chicken stock- weak/medium cereal/malty-weak	wet wall paper-medium, urine	state, wet wallpaper, musty

Results - Instrumental Analysis CF Aroma Volatiles by GC/MS

A selection of the reaction mixtures described above were further studied by GC-MS analysis of the headspace volatiles by the method described previously. These mixtures are underlined in Table 2. A sample of autolysed yeast was analysed under the same conditions for comparison of the volatiles. The results of the study are presented in detail in Table 13.

Table 13 - Volatile compounds analysed by headspace concentration

Identified Compounds	LRI	Approximate quantities (ng/0.4g of Freeze Dried Centrate, or ng/45.7g of Centrate)										
		pH5.5 140°C	ph 9 140°C	pH5.5 + C1/20 140°C	pH5.5 + C 1/5 140°C	ph 9 + C 1/20 140°C	pH 9 + C 1/5 140°C	pH5.5 + C 1/20 175°C	pH 9 + C 1/5 175°C	pH5.5 + C1/20 100°C	pH 9 + C 1/5 100°C	Yeast
Pyrazines												
pyrazine	762	-	648	-	-	1108	-	102	-	-	-	-
methylpyrazine	845	78	4607	-	37	2521	813	-	-	-	-	-
2-6-dimethylpyrazine	928	-	-	-	-	465	590	-	-	3	-	-
2,5-dimethylpyrazine	936	104	3859	-	-	382	500	-	-	17	-	4
2,3-dimethylpyrazine	940	-	966	-	-	30	-	-	-	-	-	-
ethenylpyrazine	950	44	239	-	-	-	-	-	-	157	-	-
2-ethylpyrazine	956	-	38	-	-	118	186	-	-	-	-	-
2-ethyl-6-methylpyrazine	1014	-	984	-	-	301	343	-	-	-	-	-
trimethylpyrazine	1018	256	1758	-	63	1453	1758	-	-	98	-	-
propylpyrazine	1025	-	21	-	-	-	-	-	-	-	-	-
6-methyl-2-ethenylpyrazine	1034	20	160	-	46	-	-	-	-	-	-	-
tetramethylpyrazine	1098	-	-	-	-	1009	1121	-	-	660	-	-
dimethylethylpyrazine	1099	416	120	51	-	85	111	74	-	9	-	-
dimethylethylpyrazine	1101	-	1593	-	-	-	35	-	-	-	-	-
methylpropylpyrazine	1101	26	-	-	-	-	-	-	-	-	-	-
3,5-diethyl-2-methylpyrazine	1168	-	137	-	-	-	-	18	-	-	-	4
trimethylethylpyrazine	1168	51	-	-	-	120	195	-	-	105	-	-
dimethylpropylpyrazine	1170	12	63	-	-	98	186	-	-	15	-	-
dimethylpropylpyrazine	1179	32	130	-	-	81	146	-	-	9	-	-
dimethylbutylpyrazine	1321	-	-	-	-	166	352	-	-	16	-	2
trimethylpropylpyrazine	1247	75	168	-	23	27	-	36	-	265	-	-
trimethylpropylpyrazine	1260	16	99	-	-	-	-	-	-	2	-	-
2-(2-methylpropyl)-3-(1-methylethyl)-pyrazine		-	-	-	-	-	-	-	-	-	-	-
methyldiethylpyrazine		-	-	-	-	-	-	-	-	7	-	-

2,5-dimethyl-3,6-dipropyl-pyrazine	-	-	-	-	-	-	-	-	-	5	-
Furans											
2-ethylfuran	707	-	-	29	398	32	2037	-	-	-	-
dimethylfuran	711	-	-	-	22	-	-	112	-	-	-
2,4-dimethylfuran	723	111	-	104	111	43	85	-	-	-	-
dihydro-2methyl-3(2H)-furanone	831	71	69	-	-	-	-	-	-	-	-
2-furancarboxaldehyde	857	668	-	253	109	-	-	165	-	-	-
2-furanmethanol	881	85	298	153	85	33	-	121	11	-	7
2-furanmethanethiol	926	-	-	-	249	-	-	-	406	-	-
2-acetyl-furan	924	-	-	-	-	-	-	-	-	-	-
1-(2-furanyl)-ethanone	929	84	-	629	1193	-	391	494	-	7	-
2-pentylfuran	997	113	62	?	-	35	47	-	-	-	-
1-(2-furanyl)-1-propanone	1025	21	-	-	-	-	-	-	-	-	-
Pyrans											
3,4-dihydro-2H-pyran	809	-	-	-	20	-	-	-	-	-	-
tetrahydro-6-methyl-2H-pyran- 2-one	879	-	-	-	-	-	-	-	-	8	-
tetrahydro-2H-pyran-2-oll	884	155	-	-	-	-	-	-	-	-	-
Thiophenes											
thiophene	677	-	-	-	179	-	-	-	-	-	-
methylthiophene	778	-	-	66	63	47	478	62	-	-	-
ethylthiophene	875	-	-	-	28	-	58	-	-	-	-
2,5 dimethyl thiophene	897	-	-	-	19	-	36	38	-	-	-
thiophenethiol	995	-	-	-	205	-	-	-	-	-	-
2-methyltetrahydrothiophen-3-one	1014	693	-	1366	-	-	-	-	28	-	-
2-thiophenecarboxaldehyde	1026	20	-	99	-	104	-	-	-	-	-
2-thiophenemethanethiol	1077	-	-	-	27	-	-	-	-	-	-
3-(methylthio)thiophene	1103	-	-	-	37	145	-	-	-	-	-
methylthiophene-carboxaldehyde	1143	-	-	43	5	-	-	52	97	-	-

methythiophene- carboxaldehyde	1148	-	-	526	89	134	41	329	-	-	-	-
thienothiophene	1233	-	-	35	49	-	32	16	24	-	-	-
thienothiophene	1276	-	-	-	30	-	-	-	36	-	-	-
Thiazoles												
thiazole	734	-	-	46	51	-	1165	219	870	-	-	-
methythiazole	821	-	-	-	62	72	2	71	134	-	6	-
5-methylthiazole	857	-	-	-	106	57	112	779	82	-	-	-
methythiazole	866	-	-	-	39	-	-	-	-	-	-	-
dimethylthiazole	947	-	-	-	29	88	289	-	95	-	-	-
2-ethylthiazole	959	-	-	-	22	-	42	26	78	-	-	-
2,4,5-trimethylthiazole	1016	-	-	-	1499	59	262	-	153	-	-	-
methylethylthiazole	1028	-	-	-	140	-	68	123	290	-	-	-
2-acetylthiazole	1041	-	-	-	-	647	361	-	109	-	-	-
4-methyl-5-ethenylthiazole	1043	57	67	-	-	-	59	43	-	-	20	-
2-isobutylthiazole	1065	31	-	-	-	-	-	-	-	-	-	-
methylpropylthiazole	1128	-	-	-	-	36	-	-	-	-	-	-
dimethylisopropylthiazole	1157	-	-	-	24	29	87	-	55	-	-	-
dimethylisopropylthiazole	1152	-	-	-	-	-	52	-	37	-	-	-
Aliphatic compounds												
1-propanol	616	65	-	-	101	-	-	-	-	1319	-	-
2-butanone	628	1449	-	-	2021	-	3074	-	3755	-	1164	107
ethylacetate	631	-	144	74	92	112	-	97	47	82	110	13
2,3-butanedione	637	1843	1092	2001	-	-	-	1548	-	664	-	-
3-methylbutanal	655	84	3297	1232	-	903	-	1213	752	-	-	1116
pentanal	686	2217	-	-	-	-	-	-	-	-	-	-
2-methylbutanal	691	621	703	1277	-	-	-	1242	151	-	-	582
2-methyl-1-propanol	705	1133	2412	1238	-	1614	1997	948	1791	1923	3992	17
3-methyl-2-butanone	672	-	-	-	143	-	2205	-	-	42	62	-
2-pentanal	698	-	-	48	22	-	-	34	-	-	-	-
2-pentanone	708	108	-	72	279	93	1190	85	445	34	91	-
3-pentanone	708	-	-	-	72	-	502	-	163	-	-	-
2,3-pentanedione	721	708	487	567	-	240	-	525	309	66	-	-

1-butanol	719	-	132	-	70	100	-	56	28	77	-
1-methoxy-2-propanone	723	-	-	-	-	67	-	58	-	-	-
3-hydroxy-2-butanone	728	-	-	-	1377	332	-	-	142	-	-
2-methyl-2-butenal	753	-	76	-	-	-	1694	-	-	-	-
3-methyl-3-buten-1-ol	761	511	-	1794	-	-	-	-	98	294	-
3-methylbutanol	766	682	2221	954	730	254	1375	564	355	796	94
2-butenitrile	770	-	-	-	-	-	-	-	-	-	-
heptanol	770	433	-	-	-	-	-	-	-	-	-
4-methyl-2,3-pentanedione	801	60	607	-	-	-	-	-	-	-	-
2,3-hexanedione	811	1686	137	1791	-	434	1404	-	591	-	-
hexanal	818	91	-	38	-	144	171	57	-	-	-
3-hexanone	822	36	-	-	203	-	-	244	-	58	-
2-hexanone	803	-	-	-	391	-	-	-	-	80	-
3-hydroxy-2-pentanone	831	-	-	75	59	-	125	-	-	-	-
2-heptanone	903	99	122	106	111	113	161	83	57	71	7
heptanal	910	-	-	-	-	-	-	-	-	-	138
2-hydroxy-3-hexanone	913	138	78	323	234	127	432	135	-	10	-
6-methyl-2-heptanone	967	594	722	9465	15413	1397	7347	-	-	3512	-
1,3-cyclopentanedione	973	26	24	-	-	-	-	-	-	-	-
5-methyl-2-heptanol	979	31	34	555	684	91	572	399	248	247	19
methylheptanone	984	-	-	-	149	11	-	2995	4909	110	40
octenol	988	-	-	-	-	-	-	-	-	5	23
6-methyl-5-hepten-2-one	997	-	-	114	104	-	-	30	25	20	24
2-heptenal	1036	-	-	-	-	-	-	-	-	15	81
2-nonenal	1113	-	-	-	-	-	-	70	25	7	25
2-octenal	1069	-	-	-	-	-	-	-	-	-	25
1-octanol	1078	-	-	-	-	-	-	-	-	-	-
nonanal	1116	37	-	56	-	-	48	-	-	-	-
decanal	1217	28	-	-	-	-	-	-	-	6	-
2-undecanone	1298	-	-	-	-	-	-	-	-	4	-
butanoic acid butyl ester	1380	-	-	-	-	-	-	-	-	-	3
dodecanal	1417	-	-	-	-	-	-	-	-	-	3
unsaturated ketone	1459	-	-	-	9	-	-	-	-	-	8

[illegible]

[illegible]

dimethyltrisulfide	987	-	-	402	-	2403	92	-	-	-	62	87
2-pentanethiol	841	-	-	-	-	51	-	-	-	-	-	-
1(ethylthio)-2-propanone	921	-	-	-	-	16	-	-	-	-	-	-
2,2-dithio-bis-ethanol	1280	-	-	-	-	-	-	-	-	-	-	-
Cyclic polysulfides												
3,5-dimethyl-1,2,4-trithiolane	1280	-	-	-	-	-	21	91	-	36	-	-
Dioxanes												
4-methyl-1,3-dioxane	737	321	46	256	-	-	-	-	142	-	-	-
TOTAL		18678	29633	28779	33228	16478	32517	25890	17369	11936	12195	2982

The pyrazines, thiazoles and thiophenes content was very much affected by the reaction conditions. The highest levels of pyrazines were found at pH9, as expected, since the formation of N-heterocyclic compounds in the Maillard reaction is favoured by high pH. With a few exceptions, sulphur compounds were formed only in the presence of cysteine, confirming that the content of sulphur amino acids in the freeze dried centrate was very low.

The yeast autolysate aroma volatiles were dominated by terpenes and its composition was very different from the volatiles obtained from the centrate.

Conclusion

The range of flavours was generated from the centrate, showing its potential as a flavouring ingredient or as a source of precursors for reaction product flavourings. The variables that were applied in this study were the centrate concentration, the pH, the presence of added cysteine, and the temperature/duration of the heating conditions. An addition of cysteine was necessary to generate meaty aromas which derive from sulphur-containing volatiles.

CLAIMS

1 A flavouring material for food which is an aqueous solution which comprises materials removed from a filamentous fungus by contacting it with water at an elevated temperature sufficient to reduce the nucleic acid content of the filamentous fungus in which
5 the concentration of dissolved solids is sufficiently high to render the solution stable to storage at a temperature of 20°C for a period of one month or is a solid comprising materials so removed from the filamentous fungus or is a flavouring material comprising a reaction product of such materials with a sulphur containing amino acid, hydrogen sulphide or ammonium sulphide.

10 2 A flavouring material according to Claim 1 which is an aqueous solution comprising at least 30% by weight of solids.

3 A flavouring material according to Claim 1 or 2 which comprises a reaction product of such materials with cysteine.

15 4 A process of producing a flavouring material for food which comprises the step of improving the suitability of filamentous fungus for food by subjecting it in the presence of water to a temperature sufficient to reduce its nucleic acid content thereby producing an aqueous solution and concentrating the aqueous solution.

5 5 A process as claimed in Claim 4 in which water is removed from the aqueous solution by evaporation, distillation, reverse osmosis, freeze drying, or freezing out the water as ice leaving an aqueous concentrate.

20 6 A process as claimed in Claim 4 or 5 in which the soluble components are concentrated by removing water at reduced pressure and at a temperature of 40 to 70°C.

7 A process of producing a flavouring material for food according to any one of Claims 4 to 6 in which nucleic acid is removed from a filamentous fungus in its growing state.

25 8 A process as claimed in Claim 7 in which the filamentous fungus in its growing state is washed to remove growth medium and immediately thereafter heated in the presence of water to remove nucleic acid.

9 A process as claimed in any of Claims 4 to 8 in which materials recovered from the filamentous fungus are reacted with a sulphur containing amino acid, hydrogen sulphide or
30 ammonium sulphide thereby producing savoury flavours.

10 A process as claimed in Claim 9 in which the sulphur containing amino acid is cysteine.

11 A process as claimed in Claim 9 or 10 in which the reaction is carried out in the presence of water.

35 12 A process as claimed in Claim 11 in which a solution of 1.5 to 75% of the materials removed from the fungus in water is reacted with cysteine in amounts up to 10% of cysteine by weight based on the said materials at a pH of 5.5 to 9.

13 A process as claimed in Claim 12 in which the temperature of the reaction is 110 to 140°C.

14 A food which comprises 0.1 to 15% by weight based on the solids content of a flavouring material as claimed in Claim 1, 2 or 3 or which is produced by a process as claimed in any one of Claims 4 to 13.

15 A food according to Claim 14 which is a snack, biscuits, stock, soup, stew, sauce, 5
gravy or a drink.

16 A process of flavouring food which comprises the step of improving the suitability of a filamentous fungus as food by subjecting it in the presence of water to a temperature sufficient to reduce its nucleic acid contents substantially and flavouring food with materials removed from the fungus in the said step directly or after chemical reaction.

INTERNATIONAL SEARCH REPORT

International Application No

PCT/GB 98/03716

A. CLASSIFICATION OF SUBJECT MATTER

IPC 6 A23L1/28 A23L1/227 A23L1/229

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 6 A23L

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	WO 95 23843 A (ZENECA LTD ;WARD PHILIP NIGEL (GB)) 8 September 1995 cited in the application see page 1, line 10-12; claims 1-10 ---	1-16
A	SCHINDLER F ET AL: "(In Proceedings of the International Conference on Biotechnology & Food' ((see FSTA (1991) 23 4B18)).)" FOOD BIOTECHNOLOGY, vol. 4, no. 1, 1990, pages 77-85, XP002099422 NGF Biotechnology, Hüls AG, Postfach 13 20, D-4370 Marl, Federal Republic of Germany see page 83-85 --- -/--	1-16



Further documents are listed in the continuation of box C.



Patent family members are listed in annex.

* Special categories of cited documents :

"A" document defining the general state of the art which is not considered to be of particular relevance

"E" earlier document but published on or after the international filing date

"L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)

"O" document referring to an oral disclosure, use, exhibition or other means

"P" document published prior to the international filing date but later than the priority date claimed

"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention

"X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone

"Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.

"&" document member of the same patent family

Date of the actual completion of the international search

12 April 1999

Date of mailing of the international search report

29/04/1999

Name and mailing address of the ISA

European Patent Office, P.B. 5818 Patentlaan 2
NL - 2280 HV Rijswijk
Tel. (+31-70) 340-2040, Tx. 31 651 epo nl,
Fax: (+31-70) 340-3016

Authorized officer

De Jong, E

INTERNATIONAL SEARCH REPORT

Int. l. Application No

PCT/GB 98/03716

C.(Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT

Category	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	CH 556 647 A (GREENSHIELDS INVENTIONS LTD) 13 December 1974 see column 1, line 20-60; example 1 see column 4, line 14-22 see column 8 ---	1-16
A	DATABASE WPI Section Ch, Week 8011 Derwent Publications Ltd., London, GB; Class D13, AN 80-19551C XP002098375 & JP 55 006350 B (IIZUKA C) , 15 February 1980 see abstract ---	1-16
A	PATENT ABSTRACTS OF JAPAN vol. 010, no. 362 (C-389), 4 December 1986 & JP 61 158761 A (AJINOMOTO CO INC), 18 July 1986 see abstract ---	1-16
A	NEWMARK P: "Meat substitutes. Fungal food." NATURE, UK, vol. 287, no. 5777, 1980, page 6 XP000647700 see the whole document -----	1-16

INTERNATIONAL SEARCH REPORT

Information on patent family members

International Application No

PCT/GB 98/03716

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 9523843 A	08-09-1995	AU 691576 B AU 1713395 A EP 0749471 A JP 9509580 T NZ 279872 A US 5739030 A	21-05-1998 18-09-1995 27-12-1996 30-09-1997 26-02-1998 14-04-1998
CH 556647 A	13-12-1974	NONE	