



US 20080182283A1

(19) **United States**
(12) **Patent Application Publication**
Markland et al.

(10) **Pub. No.: US 2008/0182283 A1**
(43) **Pub. Date: Jul. 31, 2008**

(54) **INHIBITORS OF HUMAN PLASMIN
DERIVED FROM THE KUNITZ DOMAINS**

(75) Inventors: **William Markland**, Southborough,
MA (US); **Robert C. Ladner**,
Ijamsville, MD (US)

Correspondence Address:
LOWRIE, LANDO & ANASTASI, LLP
ONE MAIN STREET, SUITE 1100
CAMBRIDGE, MA 02142 (US)

(73) Assignee: **DYAX CORP.**, Cambridge, MA
(US)

(21) Appl. No.: **11/928,497**

(22) Filed: **Oct. 30, 2007**

Related U.S. Application Data

(60) Continuation of application No. 11/083,742, filed on
Mar. 18, 2005, which is a continuation of application
No. 10/167,351, filed on Jun. 11, 2002, now Pat. No.
6,953,674, which is a division of application No.

09/638,770, filed on Aug. 15, 2000, now Pat. No.
6,423,498, which is a division of application No.
09/240,136, filed on Jan. 29, 1999, now Pat. No. 6,103,
499, which is a division of application No. 08/676,124,
filed on Jan. 7, 1997, now Pat. No. 6,010,880, filed as
application No. PCT/US95/00298 on Jan. 11, 1995,
which is a continuation-in-part of application No.
08/208,265, filed on Mar. 10, 1994, now abandoned,
which is a continuation-in-part of application No.
08/179,685, filed on Jan. 11, 1994, now abandoned.

Publication Classification

(51) **Int. Cl.**
C12Q 1/34 (2006.01)
(52) **U.S. Cl.** **435/18**

(57) **ABSTRACT**

This invention provides: novel proteins, which are homolo-
gous to the first Kunitz domain (K1) of lipoprotein-associated
coagulation inhibitor (LACI), and which are capable of inhib-
iting plasmin; uses of such novel proteins in therapeutic,
diagnostic, and clinical methods; and polynucleotides that
encode such novel proteins.

INHIBITORS OF HUMAN PLASMIN DERIVED FROM THE KUNITZ DOMAINS

[0001] The present application is a continuation-in-part of application Ser. No. 08/208,265 (now pending) which in turn is a continuation-in-part of Ser. No. 08/179,658, filed Jan. 11, 1994. The entirety of each of these applications is hereby incorporated by reference.

BACKGROUND OF THE INVENTION

[0002] 1. Field of the Invention

[0003] This invention relates to novel mutants of the first Kunitz domain (K_1) of the human lipoprotein-associated coagulation inhibitor LACI, which inhibit plasmin. The invention also relates to other modified Kunitz domains that inhibit plasmin and to other plasmin inhibitors.

[0004] 2. Description of the Background Art

[0005] The agent mainly responsible for fibrinolysis is plasmin, the activated form of plasminogen. Many substances can activate plasminogen, including activated Hageman factor, streptokinase, urokinase (uPA), tissue-type plasminogen activator (tPA), and plasma kallikrein (pKA). pKA is both an activator of the zymogen form of urokinase and a direct plasminogen activator.

[0006] Plasmin is undetectable in normal circulating blood, but plasminogen, the zymogen, is present at about 3 μ M. An additional, unmeasured amount of plasminogen is bound to fibrin and other components of the extracellular matrix and cell surfaces. Normal blood contains the physiological inhibitor of plasmin, α_2 -plasmin inhibitor (α_2 -PI), at about 2 μ M. Plasmin and α_2 -PI form a 1:1 complex Matrix or cell bound-plasmin is relatively inaccessible to inhibition by α_2 -PI. Thus, activation of plasmin can exceed the neutralizing capacity of α_2 -PI causing a profibrinolytic state.

[0007] Plasmin, once formed:

- [0008] i. degrades fibrin clots, sometimes prematurely;
 - [0009] ii. digests fibrinogen (the building material of clots) impairing hemostasis by causing formation of friable, easily lysed clots from the degradation products, and inhibition or platelet adhesion/aggregation by the fibrinogen degradation products;
 - [0010] iii. interacts directly with platelets to cleave glycoproteins Ib and IIb/IIIa preventing adhesion to injured endothelium in areas of high shear blood flow and impairing the aggregation response needed for platelet plug formation (ADEL86);
 - [0011] iv. proteolytically inactivates enzymes in the extrinsic coagulation pathway further promoting a prolytic state.
- [0012] Robbins (ROBB87) reviewed the plasminogen-plasmin system in detail. ROBB87 and references cited therein are hereby incorporated by reference.

Fibrinolysis and Fibrinogenolysis

[0013] Inappropriate fibrinolysis and fibrinogenolysis leading to excessive bleeding is a frequent complication of surgical procedures that require extracorporeal circulation, such as cardiopulmonary bypass, and is also encountered in thrombolytic therapy and organ transplantation, particularly liver. Other clinical conditions characterized by high incidence of bleeding diathesis include liver cirrhosis, amyloidosis, acute promyelocytic leukemia, and solid tumors. Restoration of hemostasis requires infusion of plasma and/or

plasma products, which risks immunological reaction and exposure to pathogens, e.g. hepatitis virus and HIV.

[0014] Very high blood loss can resist resolution even with massive infusion. When judged life-threatening, the hemorrhage is treated with antifibrinolytics such as ϵ -amino caproic acid (See HOOV93) (EACA), tranexamic acid, or aprotinin (NEUH89). Aprotinin is also known as Trasylol™ and as Bovine Pancreatic Trypsin Inhibitor (BPTI). Hereinafter, aprotinin will be referred to as "BPTI". EACA and tranexamic acid only prevent plasmin from binding fibrin by binding the kringles, thus leaving plasmin as a free protease in plasma. BPTI is a direct inhibitor of plasmin and is the most effective of these agents. Due to the potential for thrombotic complications, renal toxicity and, in the case of BPTI, immunogenicity, these agents are used with caution and usually reserved as a "last resort" (PUTT89). All three of the antifibrinolytic agents lack target specificity and affinity and interact with tissues and organs through uncharacterized metabolic pathways. The large doses required due to low affinity, side effects due to lack of specificity and potential for immune reaction and organ/tissue toxicity augment against use of these antifibrinolytics prophylactically to prevent bleeding or as a routine postoperative therapy to avoid or reduce transfusion therapy. Thus, there is a need for a safe antifibrinolytic. The essential attributes of such an agent are:

- [0015] i. Neutralization of relevant target fibrinolytic enzyme(s);
- [0016] ii. High affinity binding to target enzymes to minimize dose;
- [0017] iii. High specificity for target, to reduce side effects; and
- [0018] iv. High degree of similarity to human protein to minimize potential immunogenicity and organ/tissue toxicity.

All of the fibrinolytic enzymes that are candidate targets for inhibition by an efficacious antifibrinolytic are chymotrypsin-homologous serine proteases.

Excessive Bleeding

[0019] Excessive bleeding can result from deficient coagulation activity, elevated fibrinolytic activity, or a combination of the two conditions. In most bleeding diatheses one must control the activity of plasmin. The clinically beneficial effect of BPTI in reducing blood loss is thought to result from its inhibition of plasmin ($K_D \sim 0.3$ nM) or of plasma kallikrein ($K_D \sim 100$ nM) or both enzymes.

[0020] GARD93 reviews currently-used thrombolytics, saying that, although thrombolytic agents (e.g. tPA) do open blood vessels, excessive bleeding is a serious safety issue. Although tPA and streptokinase have short plasma half lives, the plasmin they activate remains in the system for a long time and, as stated, the system is potentially deficient in plasmin inhibitors. Thus, excessive activation of plasminogen can lead to a dangerous inability to clot and injurious or fatal hemorrhage. A potent, highly specific plasmin inhibitor would be useful in such cases.

[0021] BPTI is a potent plasmin inhibitor; it has been found, however, that it is sufficiently antigenic that second uses require skin testing. Furthermore, the doses of BPTI required to control bleeding are quite high and the mechanism of action is not clear. Some say that BPTI acts on plasmin while others say that it acts by inhibiting plasma kallikrein. FRAE89 reports that doses of about 840 mg of BPTI to 80 open-heart surgery patients reduced blood loss by almost half

and the mean amount transfused was decreased by 74%. Miles Inc. has recently introduced Trasylol in USA for reduction of bleeding in surgery (See Miles product brochure on Trasylol, which is hereby incorporated by reference.) LOHM93 suggests that plasmin inhibitors may be useful in controlling bleeding in surgery of the eye. SHER89 reports that BPTI may be useful in limiting bleeding in colonic surgery.

[0022] A plasmin inhibitor that is approximately as potent as BPTI or more potent but that is almost identical to a human protein domain offers similar therapeutic potential but poses less potential for antigenicity.

Angiogenesis:

[0023] Plasmin is the key enzyme in angiogenesis. OREI94 reports that a 38 kDa fragment of plasmin (lacking the catalytic domain) is a potent inhibitor of metastasis, indicating that inhibition of plasmin could be useful in blocking metastasis of tumors (FIDL94). See also ELLI92, ELLI92, OREI94 and FIDL94 and the references cited there are hereby incorporated by reference.

Plasmin

[0024] Plasmin is a serine protease derived from plasminogen. The catalytic domain of plasmin (or "CatDom") cuts peptide bonds, particularly after arginine residues and to a lesser extent after lysines and is highly homologous to trypsin, chymotrypsin, kallikrein, and many other serine proteases. Most of the specificity of plasmin derives from the kringles' binding of fibrin (LUCA83, VARA83, VARA84). On activation, the bond between ARG₅₆₁-Val₅₆₂ is cut, allowing the newly free amino terminus to form a salt bridge. The kringles remain, nevertheless, attached to the CatDom through two disulfides (COLM87, ROBB87).

[0025] BPTI has been reported to inhibit plasmin with K_D of about 300 pM (SCHN86). AUER88 reports that BPTI(R₁₅) has K_i for plasmin of about 13 nM, suggesting that R₁₅ is substantially worse than K₁₅ for plasmin binding. SCHN86 reports that BPTI in which the residues C₁₄ and C₃₈ have been converted to Alanine has K_i for plasmin of about 4.5 nM. KIDO88 reports that APP-I has K_i for plasmin of about 75 pM (7.5×10^{-11} M), the most potent inhibitor of human plasmin reported so far. DENN94a reports, however, that APP-I inhibits plasmin with $K_i = 225$ nM (2.25×10^{-7} M). Our second and third library were designed under the assumption that APP-I is a potent plasmin binder. The selection process did not select APP-I residues at most locations and the report of DENN94a explains why this happened.

[0026] With recombinant DNA techniques, it is possible to obtain a novel protein by expressing a mutated gene encoding a mutant of the native protein gene. Several strategies for picking mutations are known. In one strategy, some residues are kept constant, others are randomly mutated, and still others are mutated in a predetermined manner. This is called "variegation" and is defined in Ladner et al. U.S. Pat. No. 5,223,409, which is incorporated by reference.

[0027] DENN94a and DENN94b report selections of Kunitz domains based on APP-I for binding to the complex of Tissue Factor with Factor VII_a. They did not use LACI-K1 as parental and did not use plasmin as a target. The highest affinity binder they obtained had K_D for their target of about 2 nM. Our first-round selectants have affinity in this range, but our second round selectants are about 25-fold better than this.

[0028] Proteins taken from a particular species are assumed to be less likely to cause an immune response when injected into individuals of that species. Murine antibodies are highly antigenic in humans. "Chimeric" antibodies having human constant domains and murine variable domains are decidedly less antigenic. So called "humanized" antibodies have human constant domains and variable domains in which the CDRs are taken from murine antibodies while the framework of the variable domains are of human origin. "Humanized" antibodies are much less antigenic than are "chimeric" antibodies. In a "humanized" antibody, fifty to sixty residues of the protein are of non-human origin. The proteins of this invention comprise, in most cases, only about sixty amino acids and usually there are ten or fewer differences between the engineered protein and the parental protein. Although humans do develop antibodies even to human proteins, such as human insulin, such antibodies tend to bind weakly and the often do not prevent the injected protein from displaying its intended biological function. Using a protein from the species to be treated does not guarantee that there will be no immune response. Nevertheless, picking a protein very close in sequence to a human protein greatly reduces the risk of strong immune response in humans.

[0029] Kunitz domains are highly stable and can be produced efficiently in yeast or other host organisms. At least ten human Kunitz domains have been reported. Although APP-I was thought at one time to be a potent plasmin inhibitor, there are, actually, no human Kunitz domains that inhibit plasmin as well as does BPTI. Thus, it is a goal of the present invention to provide sequences of Kunitz domain that are both potent inhibitors of plasmin and close in sequence to human Kunitz domains.

[0030] The use of site-specific mutagenesis, whether non-random or random, to obtain mutant binding proteins of improved activity is known in the art, but success is not assured.

SUMMARY OF THE INVENTION

[0031] This invention relates to mutants of BPTI-homologous Kunitz domains that potently inhibit human plasmin. In particular, this invention relates to mutants of one domain of human LACI which are likely to be non-immunogenic to humans, and which inhibit plasmin with K_D , preferably, of about 5 nM or less, more preferably of about 300 pM or less, and most preferably about 100 pM or less. The invention also relates to the therapeutic and diagnostic use of these novel proteins.

[0032] Plasmin-inhibiting proteins are useful for the prevention or treatment of clinical conditions caused or exacerbated by plasmin, including inappropriate fibrinolysis or fibrinogenolysis, excessive bleeding associated with thrombolytics, post-operative bleeding, and inappropriate androgenesis. Plasmin-binding mutants, whether or not inhibitory, are useful for assaying plasmin in samples, in vitro, for imaging areas of plasmin activity, in vivo, and for purification of plasmin.

[0033] Preferred mutants QS4 and NS4 were selected from a library that allowed about 50 million proteins having variability at positions 13, 16, 17, 18, 19, 31, 32, 34, and 39. These proteins have an amino-acid sequence nearly identical to a human protein but inhibit plasmin with K_i of about 2 nM (i.e. about 6-fold less potent than BPTI, but 100-fold better than APP-I).

[0034] An especially preferred protein, SPI11, was selected from a library allowing variability at positions 10, 11, 13, 15, 16, 17, 18, 19, and 21 and has an affinity for plasmin which is less than 100 μ M (i.e. about 3-fold superior to BPTI in binding), and yet is much more similar in sequence to LACI, a human protein, than to the BPTI, a bovine protein. Other LACI-K1 mutants selected from this library and thought to have very high affinity for plasmin include SPI51, SPI08, and SPI23. An additional library allowing variation at positions 10, 11, 13, 15, 16, 17, 18, 19, 21, 31, 32, 34, 35, and 39 has been screened and a consensus sequence (SPIcon1) found. Variants shown to be better than QS4, and thus more preferred, include SPI51 and SPI47. Sequences that are likely to have very high affinity for plasmin yet retain an essentially human amino-acid sequence have been identified and include sequences SPI60, SPI59, SPI42, SPI55, SPI56, SPI52, SPI46, SPI49, SPI53, SPI41, and SPI57. The amino-acid sequence information that confers high affinity for the active site of plasmin can be transferred to other Kunitz domains, particularly to Kunitz domains of human origin; designs of several such proteins are disclosed.

[0035] The preferred plasmin inhibitors of the present invention fulfill one or more of the following desiderata:

[0036] 1) the K_i for plasmin is at most 20 nM, preferably not more than about 5 nM, more preferably not more than about 300 pM, and most preferably, not more than about 100 pM,

[0037] 2) the inhibitor comprises a Kunitz domain meeting the requirements shown in Table 14 with residues number by reference to BPTI,

[0038] 3) at the Kunitz domain positions 12-21 and 32-39 one of the amino-acid types listed for that position-in Table 15, and

[0039] 4) the inhibitor is more similar in amino-acid sequence to a reference sequence selected from the group SPI11, SPI15, SPI08, SPI23, SPI51, SPI47, QS4, NS4, Human LACI-K2, Human LACI-K3, Human collagen α 3 KuDom, Human TFPI-2 DOMAIN 1, Human TFPI-2 DOMAIN 2, Human TFPI-2 DOMAIN 3, HUMAN ITI-K1, Human ITI-K2, HUMAN PROTEASE NEXIN-II, Human APP-I, DPI-1.1.1, DPI-1.1.2, DPI-1.1.3, DPI-1.2.1, DPI-1.3.1, DPI-2.1, DPI-3.1.1, DPI-3.2.1, DPI-3.3.1, DPI-4.1.1, DPI-4.2.1, DPI-4.2.2, DPI-4.2.3, DPI-4.2.4, DPI-4.2.5, DPI-5.1, DPI-5.2, DPI-6.1, DPI-6.2 than is the amino acid sequence of said Kunitz domain to the sequence of BPTI.

Nomenclature

[0040] Herein, affinities are stated as $K_D(K_D(A, B)=[A]/[B]/[A-B])$. A numerically smaller K_D reflects higher affinity. For the purposes of this invention, a "plasmin inhibiting protein" is one that binds and inhibits plasmin with K_i of about 20 nM or less. "Inhibition" refers to blocking the catalytic activity of plasmin and so is measurable in vitro in assays using chromogenic or fluorogenic substrates or in assays involving macromolecules.

[0041] Amino-acid residues are discussed in three ways: full name of the amino acid, standard three-letter code, and standard single-letter code. Table use only the one-letter code. The text uses full names and three-letter code where clarity requires.

A = Ala
 C = Cys
 D = Asp
 E = Glu
 F = Phe
 G = Gly
 H = His
 I = Ile
 K = Lys
 L = Leu
 M = Met
 N = Asn
 P = Pro
 Q = Gln
 R = Arg
 S = Ser
 T = Thr
 V = Val
 W = Trp
 Y = Tyr

[0042] For the purposed of this invention, "substantially homologous" sequences are at least 51%, mmore preferably at least 80%, identical, over any specified regions. Herein, sequences that are identical are understood to be "substantially homologous". Sequences would still be "substantially homologous" if within one region of at least 20 amino acids they are sufficiently similar (51% or more) but outside the region of comparison they differed totally. An insertion of one amino acid in one sequence relative to the other counts as one mismatch. Most preferably, no more than six residues, other than at termini, are different. Preferably, the divergence in sequence, particularly in the specified regions, is in the form of "conservative modifications".

[0043] "Conservative modifications" are defined as

[0044] (a) conservative substitutions of amino acids as defined in Table 9; and

[0045] (b) single or multiple insertions or deletions of amino acids at termini, at domain boundaries, in loops, or in other segments of relatively high mobility.

[0046] Preferably, except at termini, no more than about six amino acids are inserted or deleted at any locus, and the modifications are outside regions known to contain important binding sites.

Kunitz Domains

[0047] Herein, "Kunitz domain" and "KuDom" are used interchangeably to mean a homologue of BPTI (not of the Kunitz soya-bean trypsin inhibitor). A KuDom is a domain of a protein having at least 51 amino acids (and up to about 61 amino acids) containing at least two, and preferably three, disulfides. Herein, the residues of all Kunitz domains are numbered by reference to BPTI (i.e. residues 1-58). Thus the first cysteine residue is residue 5 and the last cysteine is 55. An amino-acid sequence shall, for the purposed of this invention, be deemed a Kunitz domain if it can be aligned, with three or fewer mismatches, to the sequence shown in Table 14. An insertion or deletion of one residue shall count as one mismatch In Table 14, "x" matches any amino acid and "X" matches the types listed for that position. Disulfides bonds link at least two of: 5 to 55, 14 to 38, and 30 to 51. The number of disulfides may be reduced by one, but none of the standard cysteines shall be left unpaired. Thus, if one cysteine is changed, then a compensating cysteine is added in a suitable location or the matching cysteine is also replaced by a non-

cysteine (the latter being generally preferred). For example, *Drosophila funebris* male accessory gland protease inhibitor has no cysteine at position 5, but has a cysteine at position -1 (Gust before position 1); presumably this forms a disulfide to CYS₅₅. If Cys₁₄ and Cys₃₈ are replaced, the requirement of Gly₁₂, (Gly or Ser)₃₇, and Gly₃₆ are dropped. From zero to many residues, including additional domains (including other KuDoms), can be attached to either end of a Kunitz domain.

DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENTS

[0048] Protease inhibitors, such as Kunitz domains, function by binding into the active site of the protease so that a peptide bond (the “scissile bond”) is: 1) not cleaved, 2) cleaved very slowly, or 3) cleaved to no effect because the structure of the inhibitor prevents release or separation of the cleaved segments. In Kunitz domains, disulfide bonds act to hold the protein together even if exposed peptide bonds are cleaved. From the residue on the amino side of the scissile bond, and moving away from the bond, residues are conventionally called P1, P2, P3, etc. Residues that follow the scissile bond are called P1', P2', P3', etc. (SCHE67, SCHE68). It is generally accepted that each serine protease has sites (comprising several residues) S1, S2, etc. that receive the side groups and main-chain atoms of residues P1, P2, etc. of the substrate or inhibitor and sites S1', S2', etc. that receive the side groups and main-chain atoms of P1', P2', etc. of the substrate or inhibitor. It is the interactions between the S sites and the P side groups and main chain atoms that give the protease specificity with respect to substrates and the inhibitors specificity with respect to proteases. Because the fragment having the new amino terminus leaves the protease first, many workers designing small molecule protease inhibitors have concentrated on compounds that bind sites S1, S2, S3, etc.

[0049] LASK80 reviews protein protease inhibitors. Some inhibitors have several reactive sites on one polypeptide chains and these domains usually have different sequences, specificities, and even topologies. It is known that substituting amino acids in the P₅ to P₅' region influences the specificity of an inhibitor. Previously, attention has been focused on the P1 residue and those very close to it because these can change the specificity from one enzyme class to another. LASK80 suggests that among KuDoms, inhibitors with P1=Lys or Arg inhibit trypsin, those with P1=Tyr, Phe, Trp, Leu and Met inhibit chymotrypsin, and those with P1=Ala or Ser are likely to inhibit elastase. Among the Kazal inhibitors, LASK80 continues, inhibitors with P1=Leu or Met are strong inhibitors of elastase, and in the Bowman-Kirk family elastase is inhibited with P1=Ala, but not with P1=Leu. Such limited changes do not provide inhibitors of truly high affinity (i.e. better than 1 to 10 nM).

[0050] Kunitz domains are defined above. The 3D structure (at high resolution) of BPTI (the archetypal Kunitz domain) is known. One of the X-ray structures is deposited in the Brookhaven Protein Data Bank as “6PTI”. The 3D structure of some BPTI homologues (EIGE90, HYNE90) are known. At least seventy KuDom sequences are known. Known human homologues include three KuDoms of LACI (WUNT88, GIRA89, NOVO89), two KuDoms of Inter- α -

Trypsin Inhibitor, APP-I (KIDO88), a KuDom from collagen, and three KuDoms of TFPI-2 (SPRE94).

LACI

[0051] Lipoprotein-associated coagulation inhibitor (LACI) is a human serum phosphoglycoprotein with a molecular weight of 39 kDa (amino-acid sequence in Table 1) containing three KuDoms. We refer hereinafter to the protein as LACI and to the Kunitz domains thereof as LACI-K1 (residues 50 to 107), LACI-K2 (residues 121 to 178), and LACI-K3 (213 to 270). The cDNA sequence of LACI is reported in WUNT88. GIRA89 reports mutational studies in which the P1 residues of each of the three KuDoms were altered. LACI-K1 inhibits Factor VIIa (F.VII_a) when F.VII_a is complexed to tissue factor and LACI-K2 inhibits Factor X_a. It is not known whether LACI-K3 inhibits anything. Neither LACI nor any of the KuDoms of LACI is a potent plasmin inhibitor.

[0052] KuDoms of this invention are substantially homologous with LACI-K1, but differ in ways that confer strong plasmin inhibitory activity discussed below. Other KuDoms of this invention are homologous to other naturally-occurring KuDoms, particularly to other human KuDoms. For use in humans, the proteins of this invention are designed to be more similar in sequence to a human KuDom than to BPTI, to reduce the risk of causing an immune response.

First Library of LACI-K1 and Selectants for Binding to Plasmin

[0053] Applicants have screened a first library of LACI-K1 for mutants having high affinity for human plasmin and obtained the sequences shown in Table 2 and Table 3. These sequences may be summarized as shown in Table 16, where “preferred residues” are those appearing in at least one of the 32 variants identified as binding plasmin. The preferences at residues 13, 16, 17, 18 and 19 are strong, as shown in Table 17. Although the range of types allowed at 31 and 32 is limited, the selection indicates that an acidic group at 31 and a neutral group at 32 is preferred. At residue 17, Arg was preferred; Lys, another positively charged amino acid, was not in the library, and may be a suitable substitute for Arg. Many amino-acid types at positions 34 and 39 are consistent with high-affinity plasmin binding, but some types may hinder binding.

[0054] It should be appreciated that Applicants have not sequenced all the positive isolates of this or other libraries herein disclosed, and that some of the possible proteins may not have been present in detectable amounts.

[0055] Applicants have prepared one of the selected proteins, QS4, shown in Table 2. QS4 inhibits plasmin with a K_i of about 2 nM. Although this level of inhibition is less than that of BPTI, QS4 is a preferred molecule for use in humans because it has less potential for immunogenicity. Other proteins shown in Table 2 and Table 3 are very likely to be potent inhibitors of plasmin and are likely to pose little threat of antigenicity.

Second Library that Varies Residues 10-21

[0056] Applicants have prepared a second library of LACI-K1 derivatives shown in Table 5 and allowing variation at residues 10, 11, 13, 15, 16, 17, 18, 19, and 21. This was screened for binding to plasmin and the proteins shown in Table 6 were obtained.

[0057] “Consensus” in Table 6 is $E_{10}TGPCRARFERW_{21}$, where the seven underscored residues differ from LACI-K1. Only acidic amino acids (Glu:17 or Asp:15) were seen at position 10; Lys and Asn are not acceptable. As Glu and Asp appeared with almost equal frequency, they probably to contribute equally to binding. Acidic residues were not seen at position 11. Thr was most common (11/32) with Ser appearing often (9/32); Gly appeared 8 times. At 13, Pro was strongly preferred (24/32) with Ala second at 5/32. At 15, Arg was strongly preferred (25/32), but a few (7/32) isolates have Lys. Note that BPTI(R_{15}) is a worse plasmin inhibitor than is BPTI. At 16, Ala was preferred (22/32), but Gly did appear fairly often (10/32). At 17, Arg was most common (15/32), with Lys coming second (9/32). At residues 17 and 18, APP-I has Met and Ile. At 18, we allowed Ile or Phe. Only four isolates have Ile at 18 and none of these have Met at 17. This was surprising in view of KIDO88, but quite understandable in view of DENN94a. This collection of isolates has a broad distribution at 19: (Glu:8, Pro:7, Asp:4, Ala:3, His:3, Gly:2, Gln:2, Asn:1, Ser:1, and Arg:1), but acidic side groups are strongly preferred over basic ones. At 21, the distribution was (Trp:16, Phe:14, Leu:2, Cys:0); BPTI has Tyr at 21.

[0058] The binding of clonally pure phage that display one or another of these proteins was compared to the binding of BPTI phage (Table 6). Applicants have determined the K_i of protein SPI11 and found it to be about 88 pM which is substantially superior to BPTI.

Third Library that Varies 10-21 and 31-39

[0059] Applicants used a pool of phage of the second library (varied at residues 10, 11, 13, 15, 16, 17, 18, 19, and 21) that had been selected twice for plasmin binding as a source of DNA into which variegation was introduced at residues 31, 32, 34, 35, and 39 as shown in Table 7.

[0060] This library was screened for three rounds for binding to plasmin and the isolates shown in Table 8 were obtained. The distribution of amino-acid types is shown in Table 18 where “x” means the amino-acid type was not allowed and “*” indicates the wild-type for LACI-K1.

[0061] These sequences gave a consensus in the 10-21 and 31-40 region of $E_{10}TGPCRAKFDWR_{21} \dots E_{31}AFVYGGC \underline{GG}_{40}$ (SPIcon1 in Table 4). The ten underscored amino acids differ from LACI-K1. At eight varied positions, a second type was quite common: Asp at 10, Ala at 11, Glu at 19, Phe at 21, Thr at 31, Pro or Ser at 32, Leu or Ile at 34, and Glu at 39. At position 17, the highly potent inhibitor SPI11 has R. Thus, the sequence $D_{10}TGPCRARFDRF_{21} \dots E_{31}AFTYGGCEG_{40}$ (DPI-1.1.1 in Table 4) differs from LACI-K1 by only six residues, matches the selected sequences at the residues having strong consensus, and has preferred substitutions at positions 10, 17, 21, 34, and 39. DPI-1.1.1. is expected to have a very high affinity for plasmin and little potential for immunogenicity in humans.

[0062] Preliminary testing of proteins SPI11, BPTI, SPI23, SPI51, SPI47, QS4, SPI22, SPI54, and SPI43 for plasmin inhibitory activity placed them in the order given. SPI11 is significantly more potent than BPTI with K_i of about 88 pM. SPI23 and SPI51 are very similar in activity and only slightly less potent than BPTI. SPI47 is less potent than SPI51 but better than QS4. SPI22 is weaker than QS4. SPI54 and SPI43 are not so potent as QS4, K_i probably >4 nM.

[0063] A KuDom that is highly homologous at residues 5-55 to any one of the sequences SPI11, SPI15, SPI08, SPI23, SPI51, SPI47, QS4, and NS4, as shown in Table 4, is likely to be a potent inhibitor (K_D >5 nM) of plasmin and have a low

potential for antigenicity in humans. More preferably, to have high affinity for plasmin, a KuDom would have a sequence that is identical at residues 10-21 and 31-39 and has five or fewer differences at residues 5-9, 22-30, and 40-55 as compared to any of the sequences SPI11, SPI15, SPI08, SPI23, SPI51, SPI47, QS4, and NS4.

[0064] Using the selected sequences and the binding data of selected and natural KuDoms, we can write a recipe for a high-affinity plasmin-inhibiting KuDom that can be applied to other human KuDom parentals. First, the KuDom must meet the requirements in Table 14. The substitutions shown in Table 15 are likely to confer high-affinity plasmin inhibitory activity on any KuDom. Thus a protein that contains a sequence that is a KuDom, as shown in Table 14, and that contains at each of the position 12-21 and 32-39 an amino-acid type shown in Table 15 for that position is likely to be a potent inhibitor of human plasmin. More preferably, the protein would have an amino-acid type shown in Table 15 for all of the positions listed in Table 15. To reduce the potential for immune response, one should use one or another human KuDom as parental protein to give the sequence outside the binding region.

[0065] It is likely that a protein that comprises an amino-acid sequence that is substantially homologous to SPI11 from residue 5 through residue 55 (as shown in Table 4) and is identical to SPI11 at positions 13-19, 31, 32, 34, and 39 will inhibit human plasmin with a K_i of 5 nM or less. SPI11 differs from LACI-K1 at 7 positions. It is not clear that these substitutions are equally important in fostering plasmin binding and inhibition. There are seven molecules in which one of the substituted positions of SPI11 is changed to the residue found in LACI-K1 (i.e. “reverted”), 21 in which two of the residues are reverted, 35 in which three residues are reverted, 35 in which four are reverted, 21 in which five are reverted, and seven in which six are reverted. It is expected that those with more residues reverted will have less affinity for plasmin but also less potential for immunogenicity. A person skilled in the art can pick a protein of sufficient potency and low immunogenicity from this collection of 126. It is also possible that substitutions in SPI11 by amino acids that differ from LACI-K1 can reduce the immunogenicity without reducing the affinity for plasmin to a degree that makes the protein unsuitable for use as a drug.

Designed KuDom Plasmin Inhibitors

[0066] Hereinafter, “DPI” will mean a “Designed Plasmin Inhibitor” that are KuDoms that incorporate amino-acid sequence information from the SPI series of molecules, especially SPI11. Sequences of several DPIs and their parental proteins are given in Table 4.

[0067] Sequences DPI-1.1.1, DPI-1.1.2, DPI-1.1.3, DPI-1.1.4, DPI-1.1.5, and DPI-1.1.6 (in Table 4) differ from LACI-K1 by 6, 5, 5, 4, 3, and 2 amino acids respectively and represent a series in which affinity for plasmin may decrease slowly while similarity to a human sequence increases so as to reduce likelihood of immunogenicity. The selections from each of the libraries show that M18F is a key substitution and that either I17K or I17R is very important. Selections from the second and third library indicate that Arg is strongly preferred at 15, that an acid side group at 11 is disadvantageous to binding. The highly potent inhibitor SPI11 differs from the consensus by having R_{17} , as does BPTI. DPI-1.1.1 carries the mutations D11T, K15R, I17R, M18F, K19D, and E32A, and is likely to be highly potent as a plasmin inhibitor.

DPI-1.1.2 carries D11T, K15R, I17R, M18F, and K19D, and is likely to be highly potent. DPI-1.1.3 carries the mutations D11A, K15R, I17R, M18F, and K19D relative to LACI-K1. DPI-1.1.3 differs from DPI-1.1.2 by having A₁₁ instead of T₁₁; both proteins are likely to be very potent plasmin inhibitors. DPI-1.1.4 carries the mutations I17R, M18F, K19D, and E32A and should be quite potent. As DPI-1.1.4 has fewer of the SPI11 mutations, it may be less potent, but is also less likely to be immunogenic. DPI-1.1.5 carries the mutations I17R, M18F, and K19D. This protein is likely to be a good inhibitor and is less likely to be immunogenic. DPI-1.1.6 carries only the mutations I17R and M18F but should inhibit plasmin.

[0068] Protein DPI-1.2.1 is based on human LACI-K2 and shown in Table 4. The mutations P11T, I13P, Y17R, I18F, T19D, R32E, K34I, and L39E are likely to confer high affinity for plasmin. Some of these substitutions may not be necessary; in particular, P11T and T19D may not be necessary. Other mutations that might improve the plasmin affinity include E9A, D10E, G16A, Y21W, Y21F, R32T, K34V, and L39G.

[0069] Protein DPI-1.3.1 (Table 4) is based on human LACI-K3. The mutations R11T, L13P, N17R, E18F, N19D, R31E, P32E, K34I, and S36G are intended to confer high affinity for plasmin. Some of these substitutions may not be necessary; in particular, N19D and P32E may not be necessary. Other changes that might improve K_D include D10E, N17K, F21W and G39E.

[0070] Protein DPI-2.1 (Table 4) is based on the human collagen α 3 KuDom. The mutations E11T, T13P, D16A, F17R, I18F, L19D, A31E, R32E, and W34I are likely to confer high affinity for plasmin. Some of these substitutions may not be necessary; in particular, L19D and A31E may not be necessary. Other mutations that might improve the plasmin affinity include K9A, D10E, D16G, K20R, R32T, W34V, and G39E.

[0071] DPI-3.1.1 (Table 4) is derived from Human TFPI-2 domain 1. The exchanges Y11T, L17R, L18F, L19D, and R31E are likely to confer high affinity for plasmin. The mutation L19D may not be needed. Other mutations that might foster plasmin binding include Y21W, Y21F, Q32E, L34I, L34V, and E39G.

[0072] DPI-3.2.1 (Table 4) is derived from Human TFPI-2 domain 2. This parental domain contains insertions after residue 9 (one residue) and 42 (two residues). The mutations (V₉SVDQ₁₄ replaced by V₉ETGPC₁₄) E15R, S17K, T18F, K32T, F34V, and (H₃₉RNR₄₄ replaced by E₃₉GNR₄₄) are likely to confer affinity for plasmin. Because of the need to change the number of amino acids, DPI-3.2.1 has a higher potential for immunogenicity than do other modified human KuDoms.

[0073] DPI-3.3.1 (Table 4) is derived from human TFPI-2, domain 3. The substitutions E11T, L13P, S15R, N17R, V18F, T34I, and T36G are likely to confer high affinity for plasmin. The mutations E11T, L13P, and T34I may not be necessary. Other mutations that might foster plasmin binding include D10E, T19D, Y21W, and G39E.

[0074] DPI-4.1.1 (Table 4) is from human ITI-K1 by assertion of S10E, M15R, M17K, T18F, Q34V, and M39G. The mutations M39G and Q34V may not be necessary. Other mutations that should foster plasmin binding include: A11T, G16A, M17R, S19D, Y21W, and Y21F.

[0075] DPI-4.2.1 (Table 4) is from human ITI-K2 through the mutations V10D, R11T, F17R, I18F, and P34V. The muta-

tion P34V might not be necessary. Other mutation that should foster plasmin binding include: V10E, Q19D, L20R, W21F, P34I, and Q39E. DPI-4.2.2 is an especially preferred protein as it has only three mutations: R11T, F17R, and I18F. DPI-4.2.3 is an especially preferred protein as it has only four mutations: R11T, F17R, I18F, and L20R. DPI-4.2.4 is an especially preferred protein as it has only five mutations: R11T, F17R, I18F, L20R, and P34V. DPI-4.2.5 carries the mutations V10E, R11T, F17R, I18F, L20R, V31E, L32T, P34V, and Q39G and is highly likely to inhibit plasmin very potently. Each of the proteins DPI-4.2.1, DPI-4.2.2, DPI-4.2.3, DPI-4.2.4, and DPI-4.2.5 is very likely to be a highly potent inhibitor of plasmin.

[0076] Before DENN94a, it was thought that APP-I was a very potent plasmin inhibitor. Thus, it was surprising to select proteins from a library that was designed to allow the APP-I residues at positions 10-21 which differed strongly from APP-I. Nevertheless, APP-I can be converted into a potent plasmin inhibitor. DPI-5.1 is derived from human APP-I (also known as Protease Nexin-II) by mutations M17R and I18F and is likely to be a much better plasmin inhibitor than is APP-I itself. DPI-5.2 carries the further mutations S19D, A31E, and F34I which may foster higher affinity for plasmin.

[0077] DPI-6.1 is derived from the HKI B9 KuDom (NORR93) by the five substitutions: K11T, Q15R, T16A, M17R, and M18F. DPI-6.1 is likely to be a potent plasmin inhibitor. DPI-6.2 carries the additional mutations T19D and A34V which should foster plasmin binding.

[0078] Although BPTI is the best naturally-occurring KuDom plasmin inhibitors known, it could be improved. DPI-7.1 is derived from BPTI by the mutation I18F which is likely to increase the affinity for plasmin. DPI-7.2 carries the further mutation K15R which should increase plasmin binding. DPI-7.3 carries the added mutation R39G. DPI-7.4 carries the mutations Y10D, K15R, I18F, I19D, Q31E, and R39G and should have a very high affinity for plasmin.

Modification of Kunitz Domains

[0079] KuDoms are quite small; if this should cause a pharmacological problem, such as excessively quick elimination from circulation, two or more such domains may be joined. A preferred linker is a sequence of one or more amino acids. A preferred linker is one found between repeated domains of a human protein, especially the linkers found in human BPTI homologues, one of which has two domains (BALD85, ALBR83b) and another of which has three (WUNT88). Peptide linkers have the advantage that the entire protein may then be expressed by recombinant DNA techniques. It is also possible to use a nonpeptidyl linker, such as one of those commonly used to form immunogenic conjugates. An alternative means of increasing the serum residence of a BPTI-like KuDom is to link it to polyethyleneglycol, so called PEGylation (DAVI79).

Ways to Improve Specificity of SPI11 and other KuDom Plasmin Inhibitors:

[0080] Because we have made a large part of the surface of the KuDom SPI11 complementary to the surface of plasmin. R₁₅ is not essential for specific binding to plasmin. Many of the enzymes in the clotting and fibrinolytic pathways cut preferentially after Arg or Lys. Not having a basic residue at the P1 position may give rise to greater specificity. The variant SPI11-R15A (shown in Table 11), having an ALA at P1, is likely to be a good plasmin inhibitor and may have higher specificity for plasmin relative to other proteases than does

SPI11. The affinity of SPI11-R15A for plasmin is likely to be less than the affinity of SPI11 for plasmin, but the loss of affinity for other Arg/Lys-preferring enzymes is likely to be greater and, in many applications, specificity is more important than affinity. Other mutants that are likely to have good affinity and very high specificity include SPI11-R15G and SPI11-R15N-E32A. This approach could be applied to other high-affinity plasmin inhibitors.

Increasing the Affinity of SPI11

[0081] Variation of SPI11 as shown in Table 12 and selection of binders is likely to produce a Kunitz domain having affinity for plasmin that is higher than SPI11. This fourth library allows variegation of the 14-38 disulfide. The two segments of DNA shown are synthesized and used with primers in a PCR reaction to produce ds DNA that runs from NsiI to BstEII. The primers are identical to the 5' ends of the synthetic bits shown and of length 21 for the first and 17 for the second. As the variability is very high, we would endeavor to obtain between 10^8 and 10^9 transformants (the more the better).

Mode of Production

[0082] Proteins of this invention may be produced by any conventional technique, including

[0083] a) nonbiological synthesis by sequential coupling of components, e.g. amino acids,

[0084] b) production by recombinant DNA techniques in suitable host cells, and

[0085] c) semisynthesis, for example, by removal of undesired sequences from LACI-K1 and coupling of synthetic replacement sequences.

[0086] Proteins disclosed herein are preferably produced, recombinantly, in a suitable host, such as bacteria from the genera *Bacillus*, *Escherichia*, *Salmonella*, *Erwinia*, and yeasts from the genera *Hansenula*, *Kluyveromyces*, *Pichia*, *Rhinosporidium*, *Saccharomyces*, and *Schizosaccharomyces*, or cultured mammalian cells such as COS-1. The more preferred hosts are microorganisms of the species *Pichia pastoris*, *Bacillus subtilis*, *Bacillus brevis*, *Saccharomyces cerevisiae*, *Escherichia coli* and *Yarrowia lipolytica*. Any promoter which is functional in the host cell may be used to control gene expression.

[0087] Preferably the proteins are secreted and, most preferably, are obtained from conditioned medium. Secretion is the preferred route because proteins are more likely to fold correctly and can be produced in conditioned medium with few contaminants. Secretion is not required.

[0088] Unless there is a specific reason to include glyco-groups, we prefer proteins designed to lack N-linked glycosylation sites to reduce potential for antigenicity of glyco-groups and so that equivalent proteins can be expressed in a wide variety of organisms including: 1) *E. coli*, 2) *B. subtilis*, 3) *P. pastoris*, 4) *S. cerevisiae*, and 5) mammalian cells.

[0089] Several means exist for reducing the problem of host cells producing proteases that degrade the recombinant product; see, inter alia BANE90 and BANE91. VAND92 reports that overexpression of the *B. subtilis* signal peptidase in *E. coli* leads to increased expression of a heterologous fusion protein. ANBA88 reports that addition of PMSF (a serine proteases inhibitor) to the culture medium improved the yield of a fusion protein.

[0090] Other factors that may affect production of these and other proteins disclosed here include: 1) codon usage (optimizing codons for the host is preferred), 2) signal sequence, amino-acid sequence at intended processing sites, presence and localization of processing enzymes, deletion, mutation, or inhibition of various enzymes that might alter or degrade the engineered product and mutations that make the host more permissive in secretion (permissive secretion hosts are preferred).

[0091] Reference works on the general principles of recombinant DNA technology include Watson et al., *Molecular Biology of the Gene*, Volumes I and II, The Benjamin/Cummings Publishing Company, Inc., Menlo Park, Calif. (1987); Darnell et al., *Molecular Cell Biology*, Scientific American Books, Inc., New York, N.Y. (1986); Lewin, *Genes II*, John Wiley & Sons, New York, N.Y. (1985); Old, et al., *Principles of Gene Manipulation: An Introduction to Genetic Engineering*, 2d edition, University of California Press, Berkeley, Calif. (1981); Sambrook et al., *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Laboratory, Cold Spring Harbor, N.Y. (1989); and Ausubel et al, *Current Protocols in Molecular Biology*, Wiley Interscience, NY, (1987, 1992). These references are herein entirely incorporated by reference as are the references cited therein.

Assays for Plasmin Binding and Inhibition

[0092] Any suitable method may be used to test the compounds of this invention. Scatchard (*Ann NY Acad Sci* (1949) 51:660-669) described a classical method of measuring and analyzing binding which is applicable to protein binding. This method requires relatively pure protein and the ability to distinguish bound protein from unbound.

[0093] A second appropriate method of measuring K_D is to measure the inhibitory activity against the enzyme. If the K_D to be measured is in the 1 nM to 1 μ M range, this method requires chromogenic or fluorogenic substrates and tens of micrograms to milligrams of relatively pure inhibitor. For the proteins of this invention, having K_D in the range 5 nM to 50 pM, nanograms to micrograms of inhibitor suffice. When using this method, the competition between the inhibitor and the enzyme substrate can give a measured K_i that is higher than the true K_i . Measurement reported here are not so corrected because the correction would be very small and the any correction would reduce the K_i . Here, we use the measured K_i as a direct measure of K_D .

[0094] A third method of determining the affinity of a protein for a second material is to have the protein displayed on a genetic package, such as M13, and measure the ability of the protein to adhere to the immobilized "second material". This method is highly sensitive because the genetic packages can be amplified. We obtain at least semiquantitative values for the binding constants by use of a pH step gradient. Inhibitors of known affinity for the protease are used to establish standard profiles against which other phage-displayed inhibitors are judged. Any other suitable method of measuring protein binding may be used.

[0095] Preferably, the proteins of this invention have a K_D for plasmin of at most about 5 nM, more preferably at most about 300 pM, and most preferably 100 pM or less. Prefer-

ably, the binding is inhibitory so that K_i is the same as K_D . The K_i of QS4 for plasmin is about 2 nM. The K_i of SPI11 for plasmin is about 88 pM.

Pharmaceutical Methods and Preparations

[0096] The preferred subject of this invention is a mammal. The invention is particularly useful in the treatment of humans, but is suitable for veterinary applications too.

[0097] Herein, "protection" includes "prevention", "suppression", and "treatment". "Prevention" involves administration of drug prior to the induction of disease. "Suppression" involves administration of drug prior to the clinical appearance of disease. "Treatment" involves administration of drug after the appearance of disease.

[0098] In human and veterinary medicine, it may not be possible to distinguish between "preventing" and "suppressing" since the inductive event(s) may be unknown or latent, or the patient is not ascertained until after the occurrence of the inductive event(s). We use the term "prophylaxis" as distinct from "treatment" to encompass "preventing" and "suppressing". Herein, "protection" includes "prophylaxis". Protection need not be absolute to be useful.

[0099] Proteins of this invention may be administered, by any means, systemically or topically, to protect a subject against a disease or adverse condition. For example, administration of such a composition may be by any parenteral route, by bolus injection or by gradual perfusion. Alternatively, or concurrently, administration may be by the oral route. A suitable regimen comprises administration of an effective amount of the protein, administered as a single dose or as several doses over a period of hours, days, months, or years.

[0100] The suitable dosage of a protein of this invention may depend on the age, sex, health, and weight of the recipient, kind of concurrent treatment, if any, frequency of treatment, and the desired effect. However, the most preferred dosage can be tailored to the individual subject, as is understood and determinable by one of skill in the art, without undue experimentation by adjustment of the dose in ways known in the art.

[0101] For methods of preclinical and clinical testing of drugs, including proteins, see, e.g., Berkow et al, eds., *The Merck Manual*, 15th edition, Merck and Co., Rahway, N.J., 1987; Goodman et al., eds., *Goodman and Gilman's The Pharmacological Basis of Therapeutics*, 8th edition, Pergamon Press, Inc., Elmsford, N.Y., (1990); *Avery's Drug Treatment: Principles and Practice of Clinical Pharmacology and Therapeutics*, 3rd edition, ADIS Press, LTD., Williams and Wilkins, Baltimore, Md. (1987), Ebadi, *Pharmacology*, Little, Brown and Co., Boston, (1985), which references and references cited there are hereby incorporated by reference.

[0102] In addition to a protein here disclosed, a pharmaceutical composition may contain pharmaceutically acceptable carriers, excipients, or auxiliaries. See, e.g., Berker, supra, Goodman, supra, Avery, supra and Ebadi, supra.

In Vitro Diagnostic Methods and Reagents

[0103] Proteins of this invention may be applied in vitro to any suitable sample that might contain plasmin to measure the plasmin present. To do so, the assay must include a Signal Producing System (SPS) providing a detectable signal that depends on the amount of plasmin present. The signal may be detected visually or instrumentally. Possible signals include

production of colored, fluorescent, or luminescent products, alteration of the characteristics of absorption or emission of radiation by an assay component or product and precipitation or agglutination of a component or product.

[0104] The component of the SPS most intimately associated with the diagnostic reagent is called the "label". A label may be, e.g., a radioisotope, a fluorophore, an enzyme, a co-enzyme, an enzyme substrate, an electron-dense compound, or an agglutinable particle. A radioactive isotope can be detected by use of, for example, a γ counter or a scintillation counter or by autoradiography. Isotopes which are particularly useful are ^3H , ^{125}I , ^{131}I , ^{35}S , ^{14}C , and, preferably, ^{125}I . It is also possible to label a compound with a fluorescent compound. When the fluorescently labeled compound is exposed to light of the proper wave length, its presence can be detected. Among the most commonly used fluorescent labeling compounds are fluorescein isothiocyanate, rhodamine, phycoerythrin, phycocyanin, allophycocyanin, o-phthaldehyde, and fluorescamine. Alternatively, fluorescence-emitting metals, such as ^{125}Eu or other lanthanide, may be attached to the binding protein using such metal chelating groups as diethylenetriaminepentaacetic acid or ethylenediamine-tetraacetic acid. The proteins also can be detectably labeled by coupling to a chemiluminescent compound, such as luminol, isolumino, therromatic acridinium ester, imidazole, acridinium salt, and oxalate ester. Likewise, a bioluminescent compound, such as luciferin, luciferase and aequorin, may be used to label the binding protein. The presence of a bioluminescent protein is determined by detecting the presence of luminescence. Enzyme labels, such as horseradish peroxidase and alkaline phosphatase, are preferred.

[0105] There are two basic types of assays: heterogeneous and homogeneous. In heterogeneous assays, binding of the affinity molecule to analyte does not affect the label; thus, to determine the amount of analyte, bound label must be separated from free label. In homogeneous assays, the interaction does affect the activity of the label, and analyte can be measured without separation.

[0106] In general, a plasmin-binding protein (PBP) may be used diagnostically in the same way that an antiplasmin antibody is used. Thus, depending on the assay format, it may be used to assay plasmin, or, by competitive inhibition, other substances which bind plasmin.

[0107] The sample will normally be a biological fluid, such as blood, urine, lymph, semen, milk, or cerebrospinal fluid, or a derivative thereof, or a biological tissue, e.g., a tissue section or homogenate. The sample could be anything. If the sample is a biological fluid or tissue, it may be taken from a human or other mammal, vertebrate or animal, or from a plant. The preferred sample is blood, or a fraction or derivative thereof.

[0108] In one embodiment, the plasmin-binding protein (PBP) is immobilized, and plasmin in the sample is allowed to compete with a known quantity of a labeled or specifically labelable plasmin analogue. The "plasmin analogue" is a molecule capable of competing with plasmin for binding to the PBP, which includes plasmin itself. It may be labeled already, or it may be labeled subsequently by specifically binding the label to a moiety differentiating the plasmin analogue from plasmin. The phases are separated, and the labeled plasmin analogue in one phase is quantified.

[0109] In a "sandwich assay", both an insolubilized plasmin-binding agent (PBA), and a labeled PBA are employed. The plasmin analyte is captured by the insolubilized PBA and is tagged by the labeled PBA, forming a tertiary complex. The

reagents may be added to the sample in any order. The PBAs may be the same or different, and only one PBA need be a PBP according to this invention (the other may be, e.g., an antibody). The amount of labeled PBA in the tertiary complex is directly proportional to the amount of plasmin in the sample.

[0110] The two embodiments described above are both heterogeneous assays. A homogeneous assay requires only that the label be affected by the binding of the PBP to plasmin. The plasmin analyte may act as its own label if a plasmin inhibitor is used as a diagnostic reagent.

[0111] A label may be conjugated, directly or indirectly (e.g., through a labeled anti-PBP antibody), covalently (e.g., with SPDP) or noncovalently, to the plasmin-binding protein, to produce a diagnostic reagent. Similarly, the plasmin binding protein may be conjugated to a solid phase support to form a solid phase ("capture") diagnostic reagent. Suitable supports include glass, polystyrene, polypropylene, polyethylene, dextran, nylon, amylases, and magnetite. The carrier can be soluble to some extent or insoluble for the purposes of this invention. The support material may have any structure so long as the coupled molecule is capable of binding plasmin.

In Vivo Diagnostic Uses

[0112] A Kunitz domain that binds very tightly to plasmin can be used for in vivo imaging. Diagnostic imaging of disease foci was considered one of the largest commercial opportunities for monoclonal antibodies, but this opportunity has not been achieved. Despite considerable effort, only two monoclonal antibody-based imaging agents have been approved. The disappointing results obtained with monoclonal antibodies is due in large measure to:

[0113] i) Inadequate affinity and/or specificity;

[0114] ii) Poor penetration to target sites;

[0115] iii) Slow clearance from nontarget sites;

[0116] iv) Immunogenicity (most are murine); and

[0117] v) High production cost and poor stability.

[0118] These limitations have led most in the diagnostic imaging field to begin to develop peptide-based imaging agents. While potentially solving the problems of poor penetration and slow clearance, peptide-based imaging agents are unlikely to possess adequate affinity, specificity and in vivo stability to be useful in most applications.

[0119] Engineered proteins are uniquely suited to the requirements for an imaging agent. In particular the extraordinary affinity and specificity that is obtainable by engineering small, stable, human-origin protein domains having known in vivo clearance rates and mechanisms combine to provide earlier, more reliable results, less toxicity/side effects, lower production and storage cost, and greater convenience of label preparation. Indeed, it should be possible to achieve the goal of realtime imaging with engineered protein imaging agents. Plasmin-binding proteins, e.g. SPI11, may be useful for localizing sites of internal hemorrhage.

[0120] Radio-labelled binding protein may be administered to the human or animal subject. Administration is typically by injection, e.g., intravenous or arterial or other means of administration in a quantity sufficient to permit subsequent dynamic and/or static imaging using suitable radio-detecting devices. The dosage is the smallest amount capable of providing a diagnostically effective image, and may be determined by means conventional in the art, using known radio-imaging agents as guides.

[0121] Typically, the imaging is carried out on the whole body of the subject, or on that portion of the body or organ relevant to the condition or disease under study. The radio-labelled binding protein has accumulated. The amount of radio-labelled binding protein accumulated at a given point in time in relevant target organs can then be quantified.

[0122] A particularly suitable radio-detecting device is a scintillation cameras, such as a γ camera. The detection device in the camera senses and records (and optional digitizes) the radioactive decay. Digitized information can be analyzed in any suitable way, many of which are known in the art. For example, a time-activity analysis can illustrate uptake through clearance of the radio-labelled binding protein by the target organs with time.

[0123] Various factors are taken into consideration in picking an appropriate radioisotope. The isotope is picked: to allow good quality resolution upon imaging, to be safe for diagnostic use in humans and animals, and, preferably, to have a short half-life so as to decrease the amount of radiation received by the body. The radioisotope used should preferably be pharmacologically inert, and the quantities administered should not have substantial physiological effect. The binding protein may be radio-labelled with different isotopes of iodine, for example ^{123}I , ^{125}I , or ^{131}I (see, for example, U.S. Pat. No. 4,609,725). The amount of labeling must be suitably monitored.

[0124] In applications to human subjects, it may be desirable to use radioisotopes other than ^{125}I for labelling to decrease the total dosimetry exposure of the body and to optimize the detectability of the labelled molecule. Considering ready clinical availability for use in humans, preferred radio-labels include: ^{99m}Tc , ^{67}Ga , ^{68}Ga , ^{90}Y , ^{111}In , ^{113m}In , ^{123}I , ^{186}Re , ^{188}Re or ^{211}At . Radio-labelled protein may be prepared by various methods. These include radio-halogenation by the chloramine-T or lactoperoxidase method and subsequent purification by high pressure liquid chromatography, for example, see Gutkowska et al in "Endocrinology and Metabolism Clinics of America: (1987) 16 (1):183. Other methods of radio-labelling can be used, such as IODO-BEADSTM.

[0125] A radio-labelled protein may be administered by any means that enables the active agent to reach the agent's site of action in a mammal. Because proteins are subject to digestion when administered orally, parenteral administration, i.e., intravenous subcutaneous, intramuscular, would ordinarily be used to optimize absorption.

Other Uses

[0126] The plasmin-binding proteins of this invention may also be used to purify plasmin from a fluid, e.g., blood. For this purpose, the PBP is preferably immobilized on an insoluble support. Such supports include those already mentioned as useful in preparing solid phase diagnostic reagents.

[0127] Proteins can be used as molecular weight markers for reference in the separation or purification of proteins. Proteins may need to be denatured to serve as molecular weight markers. A second general utility for proteins is the use of hydrolyzed protein as a nutrient source. Proteins may also be used to increase the viscosity of a solution.

[0128] The protein of his invention may be used for any of the foregoing purposes, as well as for therapeutic and diagnostic purposes as discussed further earlier in this specification.

Preparation of Peptides

[0129] Chemical polypeptide synthesis is a rapidly evolving area in the art, and methods of solid phase polypeptide synthesis are well-described in the following references, hereby entirely incorporated by reference: (Merrifield, *J Amer Chem Soc* 85:2149-2154 (1963); Merrifield, *Science* 232:341-347 (1986); Wade et al., *Biopolymers* 25:S21-S37 (1986); Fields, *Int J Polypeptide Prot Res* 35:161 (1990); MilliGen Report Nos. 2 and 2a, Millipore Corporation, Bedford, Mass., 1987) Ausubel et al, supra, and Sambrook et al, supra. Tan and Kaiser (*Biochemistry*, 1977, 16:1531-41) synthesized BPTI and a homologue eighteen years ago.

[0130] As is known in the art, such methods involve blocking or protecting reactive functional groups, such as free amino, carboxyl and thio groups. After polypeptide bond formation, the protective groups are removed. Thus, the addition of each amino acid residue requires several reaction steps for protecting and deprotecting. Current methods utilize solid phase synthesis, wherein the C-terminal amino acid is covalently linked to an insoluble resin particles that can be filtered. Reactants are removed by washing the resin particles with appropriate solvents using an automated machine. Various methods, including the "tBoc" method and the "Fmoc" method are well known in the art. See, inter alia, Atherton et al., *J Chem Soc Perkin Trails* 1:538-546 (1981) and Sheppard et al., *Int J Polypeptide Prot Res* 20:451-454 (1982).

EXAMPLES

Example 1

Construction of LACI (K1) Library

[0131] A synthetic oligonucleotide duplex having NsiI- and MluI-compatible ends was cloned into a parental vector (LACI-K1::IU) previously cleaved with the above two enzymes. The resultant ligated material was transfected by electroporation into XLIMR (F⁻) *E. coli* strain and plated on ampicillin (Ap) plates to obtain phage-generating Ap^R colonies. The variegation scheme for Phase I focuses on the P1 region, and affected residues 13, 16, 17, 18 and 19. It allowed for 6.6×10⁵ different DNA sequences (3.1×10⁵ different protein sequences). The library obtained consisted of 1.4×10⁶ independent cfu's which is approximately a two fold representation of the whole library. The phage stock generated from this plating gave a total titer of 1.4×10¹³ pfu's in about 3.9 ml, with each independent done being represented, on average, 1×10⁷ in total and 2.6×10⁶ times per ml of phage stock.

[0132] To allow for variegation of residues 31, 32, 34 and 39 (phase II), synthetic oligonucleotide duplexes with MluI- and BstEII-compatible ends were cloned into previously cleaved R_γ-DNA derived from one of the following

[0133] i) the parental construction,

[0134] ii) the phase I library, or

[0135] iii) display phage selected from the first phase binding to a given target.

[0136] The variegation scheme for phase II allows for 4096 different DNA sequences (1600 different protein sequences) due to alterations at residues 31, 32, 34 and 39. The final phase

II variegation is dependent upon the level of variegation remaining following the three rounds of binding and elution with a given target in phase I.

[0137] The combined possible variegation for both phases equals 2.7×10⁸ different DNA sequences or 5.0×10⁷ different protein sequences. When previously selected display phage are used as the origin of R_γ-DNA for the phase II variegation, the final level of variegation is probably in the range of 10⁵ to 10⁶.

Example 2

Screening of LACI-K1 Library for Binding to Plasmin

[0138] The scheme for selecting LACI-K1 variants that bind plasmin involves incubation of the phage-display library with plasmin-beads (Calbiochem, San Diego, Calif.; catalogue no. 527802) in a buffer (PBS containing 1 mg/ml BSA) before washing away unbound and poorly retained display-phage variants with PBS containing 0.1% Tween 20. The more strongly bound display-phage are eluted with a low pH elution buffer, typically citrate buffer (pH 2.0) containing 1 mg/ml BSA, which is immediately neutralized with Tris buffer to pH 7.5. This process constitutes a single round of selection.

[0139] The neutralized eluted display-phage can be either used:

[0140] i) to inoculate an F⁺ strain of *E. coli* to generate a new display-phage stock, to be used for subsequent rounds of selection (so-called conventional screening), or

[0141] ii) be used directly for another immediate round of selection with the protease beads (so-called quick screening).

[0142] Typically, three rounds of either method, or a combination of the two, are performed to give rise to the final selected display-phage from which a representative number are sequenced and analyzed for binding properties either as pools of display-phage or as individual clones.

[0143] For the LACI-K1 library, two phases of selection were performed, each consisting of three rounds of binding and elution. Phase I selection used the phase I library (variegated residues 13, 16, 17, 18, and 19) which went through three rounds of binding and elution against plasmin giving rise to a subpopulation of clones. The R_γ-DNA derived from this selected subpopulation was used to generate the Phase II library (addition of variegated residues 31, 32, 34 and 39). About 5.6×10⁷ independent transformants were obtained. The phase II libraries underwent three further rounds of binding and elution with the same target protease giving rise to the final selectants.

[0144] Following two phases of selection against plasmin-agarose beads a representative number (16) of final selection display-phage were sequenced. Table 2 shows the sequences of the selected LACI-K1 domains with the amino acids selected at the variegated positions in upper case. Note the absolute selection of residues P₁₃, A₁₆, R₁₇, F₁₁, and E₁₉. There is very strong selection for E at 31 and Q at 32. There is no consensus at 34; the observed amino acids are {T₃, Y₂, H₂, D, R, A, V₂, I₃, and L}. The amino acids having side groups that branch at C₀ (T, I, and V) are multiply represented and are preferred. At position 39, there is no strong consensus (G₆, D₃, Q₂, A₂, R, F, E), but G, D, Q, and A seem to be preferred (in that order).

[0145] A separate screening of the LACI-K1 library against plasmin gave a very similar consensus from 16 sequenced selected display-phage. These sequences are shown in Table 3 (selected residues in upper case). These sequences depart from those of Table 2 in that E here predominates at position 19. There is a consensus at 34 (T₃, V₃, S₃, I₂, L, A, F) of T, V, or S. Combining the two sets, there is a preference for (in order of preference) T, V, I, S, A, H, Y, and L, with F, D, and R being allowed.

Expression, Purification and Kinetic Analysis.

[0146] The three isolates QS4, ARFK#1, and ARFK#2 were recloned into a yeast expression vector. The yeast expression vector is derived from pMFalpha8 (KURJ82 and MIYA85). The LACI variant genes were fused to part of the matα 1 gene, generating a hybrid gene consisting of the matα 1 promoter-signal peptide, and leader sequence-fused to the LACI variant. The cloning site is shown in Table 24. Note that the correctly processed LACI-K1 variant protein should be as detailed in Table 2 and Table 3 with the addition of residues glu-ala-ala-glu to the N-terminal met (residue 1 in Table 2 and Table 3). Expression in *S. cerevisiae* gave a yield of about 500 μg of protease inhibitor per liter of medium. Yeast-expressed LACI (kunitz domain 1), BPTI and LACI variants: QS4, ARFK#1 and ARFK#2 were purified by affinity chromatography using trypsin-agarose beads.

[0147] The most preferred production host is *Pichia pastoris* utilizing the alcohol oxidase system Others have produced a number of proteins in the yeast *Pichia pastoris*. For example, Vedvick et al. (VEDV91) and Wagner et al. (WAGN92) produced aprotinin from the alcohol oxidase promoter with induction by methanol as a secreted protein in the culture medium at ≈1 mg/ml. Gregg et al. (GREG93) have reviewed production of a number of proteins in *P. pastoris*. Table 1 of GREG93 shows proteins that have been produced in *P. pastoris* and the yields.

Kinetic Data

[0148] Inhibition of hydrolysis of succinyl-Ala-Phe-Lys-(F₃Ac)AMC (a methyl coumarin) (Sigma Chemical, St. Louis, Mo.) by plasmin at 2.5×10⁻¹ M with varying amount of inhibitor were fit to the standard form for a tight-binding substrate by least-squares. Preliminary kinetic analysis of the two ARFK variants demonstrated very similar inhibitory activity to that of the QS4 variant.) These measurements were carried out with physiological amounts of salt (150 mM) so that the affinities are relevant to the action of the proteins in blood.

[0149] Table 23 shows that QS4 is a highly specific inhibitor of human plasmin. Phage that display the LACI-K1

derivative QS4 bind to plasmin beads at least 50-times more than it binds to other protease targets.

New Library for Plasmin:

[0150] A new library of LACI-K1 domains, displayed on M13 gIIIp and containing the diversity shown in Table 5 was made and screened for plasmin binding. Table 6 shows the sequences selected and the consensus. We characterized the binding of the selected proteins by comparing the binding of clonally pure phage to BPTI display phage. Isolates 11, 15, 08, 23, and 22 were superior to BPTI phage. We produced soluble SPI11(Selected Plasmin Inhibitor #11) and tested its inhibitory activity, obtaining a K_i of 88 pM which is at least two-fold better than BPTI. Thus, we believe that the selectants SPI15, SPI08, and SPI22 are far superior to BPTI and that SPI23 is likely to be about as potent as BPTI. All of the listed proteins are much closer to a human protein amino-acid sequence than is BPTI and so have less potential for immunogenicity.

[0151] All references, including those to U.S. and foreign patents or patent applications, and to nonpatent disclosures, are hereby incorporated by reference in their entirety.

TABLE 1

Sequence of whole LACI: (SEQ ID NO. 1)			
	5	5	5
1	MIYTHKKVHA	LWASVCLLLN	LAPAPLNAds
	5		5
	edehehtiit		
	dtelpplkIM		
51	HSFCAFKADD	GPCKAIMKRF	FFNIFTRQCE
	QNRFSLEEC		EFIYGGCEGN
101	KKMCTRDnan	riiktntlqqe	kpdcfleeed
	yfyngtkqC		pgicrgvitr
151	erfkyqgClq	nmnnfetlee	CkniCedgpn
	lnavnsltp		gfvdnvqtq
201	qstkvpslfe	fhgpswCltp	adrglCrane
	krpfkysgC		nrfyynsvig
251	ggnennftsk	qeClraCkkg	fiqriskggl
	rvkiayeeif		iktkrkrkkq
301	vknm		

The signal sequence (1-28) is uppercase and underscored
 LACI-K1 is uppercase
 LACI-K2 is underscored
 LACI-K3 is bold

TABLE 2

Sequence of LACI-K1 and derivatives that bind human plasmin					
	1	2	3	4	5
	1234567890123456789012345678901234567890123456789012345678				
LACI-K1	mhsfcaf	kaddg	pckaimkr	ffniftrqce	efiyggcegnqnr
	fesleeck	kkmctrd			SEQ ID NO. 2
QS1	mhsfcaf	kaddg	PckARF	ErffniftrqcEQf	TyggcRgnqnr
	fesleeck	kkmctrd			SEQ ID NO. 3
QS4	mhsfcaf	kaddg	PckARF	ErffniftrqcEQf	YggcDgnqnr
	fesleeck	kkmctrd			SEQ ID NO. 4

TABLE 2-continued

<u>Sequence of LACI-K1 and derivatives that bind human plasmin</u>						
	1	2	3	4	5	
	1234567890123456789012345678901234567890123456789012345678					
QS7	mhsfcafkaddgPckARFERfffniftrqcEQfHyggcDgnqnrfesleeckkmctrd					SEQ ID NO. 5
QS8	mhsfcafkaddgPckARFERfffniftrqcEQfDyggcAgnqnrfesleeckkmctrd					SEQ ID NO. 6
QS9	mhsfcafkaddgPckARFERfffniftrqcEQfRyggcDgnqnrfesleeckkmctrd					SEQ ID NO. 7
QS13	mhsfcafkaddgPckARFERfffniftrqcQQfYyggcQgnqnrfesleeckkmctrd					SEQ ID NO. 8
QS15	mhsfcafkaddgPckARFERfffniftrqcEEfAyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 9
NS2	mhsfcafkaddgPckARFERfffniftrqcQQfVyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 10
NS4	mhsfcafkaddgPckARFERfffniftrqcEQfTyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 11
NS6	mhsfcafkaddgPckARFERfffniftrqcEEfTyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 12
NS9	mhsfcafkaddgPckARFERfffniftrqcEQfIyggcQgnqnrfesleeckkmctrd					SEQ ID NO. 13
NS11	mhsfcafkaddgPckARFERfffniftrqcEQfIyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 14
NS12	mhsfcafkaddgPckARFERfffniftrqcEQfIyggcFgnqnrfesleeckkmctrd					SEQ ID NO. 15
NS14	mhsfcafkaddgPckARFERfffniftrqcQQfHyggcEgnqnrfesleeckkmctrd					SEQ ID NO. 16
NS15	mhsfcafkaddgPckARFERfffniftrqcEQfVyggcAgnqnrfesleeckkmctrd					SEQ ID NO. 17
NS16	mhsfcafkaddgPckARFERfffniftrqcEQfLyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 18
ARFRCON	mhsfcafkaddgPckARFERfffniftrqcEQfiyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 19

TABLE 3

<u>Sequence of LACI-K1 and derivatives that bind human plasmin</u>						
	1	2	3	4	5	
	1234567890123456789012345678901234567890123456789012345678					
LACI-K1	mhsfcafkaddgpckaimkrfffniftrqceefiyggcegnqnrfesleeckkmctrd					SEQ ID NO. 2
ARFK#1	mhsfcafkaddgPckARFERfffniftrqcEQfVyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 20
ARFK#2	mhsfcafkaddgPckARFERfffniftrqcEEfVyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 21
ARFK#3	mhsfcafkaddgLckGRFQRfffniftrqcEEfIyggcEgnqnrfesleeckkmctrd					SEQ ID NO. 22
ARFK#4	mhsfcafkaddgPckARFERfffniftrqcEQfTyggcMgnqnrfesleeckkmctrd					SEQ ID NO. 23
ARFK#5	mhsfcafkaddgPckARFERfffniftrqcEQfSyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 24
ARFK#6	mhsfcafkaddgPckARFERfffniftrqcEEfLyggcLgnqnrfesleeckkmctrd					SEQ ID NO. 25
ARFK#7	mhsfcafkaddgPckARFERfffniftrqcEQfSyggcQgnqnrfesleeckkmctrd					SEQ ID NO. 26
ARFK#8	mhsfcafkaddgPckARFERfffniftrqcEQfAyggcAgnqnrfesleeckkmctrd					SEQ ID NO. 27
ARFK#9	mhsfcafkaddgPckARFERfffniftrqcEQfIyggcVgnqnrfesleeckkmctrd					SEQ ID NO. 28
ARFK#10	mhsfcafkaddgPckARFERfffniftrqcEEfSyggcKgnqnrfesleeckkmctrd					SEQ ID NO. 29
ARFK#11	mhsfcafkaddgPckARFERfffniftrqcEEfVyggcKgnqnrfesleeckkmctrd					SEQ ID NO. 30
ARFK#12	mhsfcafkaddgPckASFERfffniftrqcEQfTyggcNgnqnrfesleeckkmctrd					SEQ ID NO. 31
ARFK#13	mhsfcafkaddgPckASFERfffniftrqcEEfTyggcLgnqnrfesleeckkmctrd					SEQ ID NO. 32
ARFK#14	mhsfcafkaddgPckARFERfffniftrqcEQfVyggcHgnqnrfesleeckkmctrd					SEQ ID NO. 33
ARFK#15	mhsfcafkaddgPckARFERfffniftrqcEQfTyggcGgnqnrfesleeckkmctrd					SEQ ID NO. 34

TABLE 3-continued

<u>Sequence of LACI-K1 and derivatives that bind human plasmin</u>	
1	2
3	4
5	
123456789012345678901234567890123456789012345678	
ARFK#16	mhsfcfakaddgPckARFERfffniftrqcEQfTyggcMgnqnrfsleeckkmctrd SEQ ID NO. 35
ARFKC01	mhsfcfakaddgPckARFERfffniftrqcEQfTyggcGgnqnrfsleeckkmctrd SEQ ID NO. 36
ARFKC02	mhsfcfakaddgPckARFERfffniftrqcEQfVyggcGgnqnrfsleeckkmctrd SEQ ID NO. 37

TABLE 4

<u>Kunitz domains, some of which inhibit plasmin</u>			
Protein identifier	Amino-acid Sequence	Affinity K _D	SEQ ID NO.
QS4	mhsfcfakaddgPckARFERfffniftrqcEQfTyggcDgnqnrfsleeckkmctrd	2 nM	SEQ ID NO. 4
NS4	mhsfcfakaddgPckARFERfffniftrqcEQfTyggcGgnqnrfsleeckkmctrd	(B)	SEQ ID NO. 11
BPTI	RPDFCLEPPYTPGPKARIIRYFYNAKAGLCQTFVYGGCRAKRNNFKSAEDCMRTCGGA	.3 nM	SEQ ID NO. 38
Human APP-I	VREVCSEQAETGPCRAMISRWFYFDVTEGKCAPFFYGGCGGNRNNFDTEEYCMAVCGSA	75 pM (KIDO88) 225 nM (DENN94a)	SEQ ID NO. 39
SpI11	mhsfcfakaETgPcRARFDrWffniftrqceefiyggcegnqnrfsleeckkmctrd	88 pM	SEQ ID NO. 40
SpI15	mhsfcfakaESgPcRARFDrWffniftrqceefiyggcegnqnrfsleeckkmctrd	(A)	SEQ ID NO. 41
SPI08	mhsfcfakaDGgPcRARFERfffniftrqceefiyggcegnqnrfsleeckkmctrd	(A)	SEQ ID NO. 42
SPI23	mhsfcfakaEGgPcRAKfQrWffniftrqceefiyggcegnqnrfsleeckkmctrd	~.5 nM	SEQ ID NO. 43
SPI22	mhsfcfakaDGgPcKGFPrFffniftrqceefiyggcegnqnrfsleeckkmctrd	>2 nM	SEQ ID NO. 44
SPIcon1	mhsfcfakaETgPcRAkFDrWffniftrqcEAFVyggcGgnqnrfsleeckkmctrd		SEQ ID NO. 45
SPI60	mhsfcfakaETgPcRAkFDrWffniftrqcEPfVYggcEgnqnrfsleeckkmctrd	(B)	SEQ ID NO. 46
SPI59	mhsfcfakaETgPcRAkFDrWffniftrqcNTfVYggcGgnqnrfsleeckkmctrd		SEQ ID NO. 47
SPI42	mhsfcfakaETgPcRGkFDrWffniftrqcQGfVYggcGgnqnrfsleeckkmctrd		SEQ ID NO. 48
SPI55	mhsfcfakaEVgPcRAkFDrWffniftrqcHLfTYggcGgnqnrfsleeckkmctrd		SEQ ID NO. 49
SPI56	mhsfcfakaETgPcRGkFDrWffniftrqcAQfVYggcEgnqnrfsleeckkmctrd		SEQ ID NO. 50
SPI43	mhsfcfakaETgPcRGkFDrWffniftrqcESfHYggcKgnqnrfsleeckkmctrd	>-4 nM	SEQ ID NO. 51
SPI52	mhsfcfakaDAGpCRAkFERfffniftrqcEAFLYggcGgnqnrfsleeckkmctrd		SEQ ID NO. 52
SPI46	mhsfcfakaDVgPcRAkFERfffniftrqcEAFLYggcEgnqnrfsleeckkmctrd		SEQ ID NO. 53
SPI51	mhsfcfakaDAGpCRAkFERfffniftrqcTAFYggcGgnqnrfsleeckkmctrd	~.5 nM	SEQ ID NO. 54
SPI54	mhsfcfakaDSgPcRARFDrWffniftrqcTRfPYggcGgnqnrfsleeckkmctrd	>-4 nM	SEQ ID NO. 55
SPI49	mhsfcfakaETgPcRAkIPrLffniftrqcEPfIWggcGgnqnrfsleeckkmctrd		SEQ ID NO. 56
SPI47	mhsfcfakaDAGpCRAkFERfffniftrqcEEfIYggcEgnqnrfsleeckkmctrd	-.8 nM	SEQ ID NO. 57
SPI53	mhsfcfakaETgPcKGSFDrWffniftrqcNVfRYggcRgnqnrfsleeckkmctrd		SEQ ID NO. 58
SPI41	mhsfcfakaDAGpCRAkFERfffniftrqcDTfLYggcEgnqnrfsleeckkmctrd	(AB)	SEQ ID NO. 59
SPI57	mhsfcfakaDSgPcKGRFGrLffniftrqcTAFDWggcGgnqnrfsleeckkmctrd		SEQ ID NO. 60

TABLE 4-continued

<u>Kunitz domains, some of which inhibit plasmin</u>			
Protein identifier	Amino-acid Sequence 1111111111222222222233333333334444444444A5555555555	Affinity K _D	SEQ ID NO.
DPI-1.1.1	mhsfcafkadTgpcRaRFDrfffniftrqceAfiyggcegnqnrfsleeckkmctrd	(A)	SEQ ID NO. 61
DPI-1.1.2	mhsfcafkadTgpcRaRFDrfffniftrqceefiyggcegnqnrfsleeckkmctrd	(A)	SEQ ID NO. 62
DPI-1.1.3	mhsfcafkadAgpcRaRFDrfffniftrqceefiyggcegnqnrfsleeckkmctrd	(A)	SEQ ID NO. 63
DPI-1.1.4	mhsfcafkadTgpcKaRFDrfffniftrqceAfiyggcegnqnrfsleeckkmctrd	(AB)	SEQ ID NO. 127
DPI-1.1.5	mhsfcafkaddgpckaRFDrfffniftrqceefiyggcegnqnrfsleeckkmctrd	(B)	SEQ ID NO. 128
DPI-1.1.6	mhsfcafkaddgpckaRfkrfffniftrqceefiyggcegnqnrfsleeckkmctrd	(C)	SEQ ID NO. 129
Human LACI-K2	KPDFCFLEEDPGICRGYITRYFYNNQTKQCFKYGGLGNNMNFETLEECKNICEDG		SEQ ID NO. 64
DPI-1.2.1	kpdfcfleedTgPcrgRFDryfynnqtkqceTfIyggcEgnmnnfetleecknicedg		SEQ ID NO. 65
Human LACI-K3	GPSWCLTPADRGLCRANENRFYNSVIGKCRPFKYSGCCGNNENFTSKQECLRACKKG		SEQ ID NO. 66
DPI-1.3.1	gpswcltpadTgPcraRFDryfynsvigkcEpfiyGgcegnennftskqeclrackkg		SEQ ID NO. 67
Human collagen α3 KuDom	ETDICKLPKDEGTCRDFILKWYDPNTKSCARFWYGGCGGNNKFGSQKECEKVCAPV		SEQ ID NO. 68
DPI-2.1	etdicklpkdTgPcraRFDkwydpntkscEPfVyggegnenkfsgsqkecekvcapv		SEQ ID NO. 69
Human TFPI-2 DOMAIN 1	NAEICLLPLDYGPCRALLLRYYDRYTQSCRQFLYGGCEGNANNFYTWEACDDACWRI		SEQ ID NO. 70
DPI-3.1.1	naeicllpldTgpcraRFDryyydrytqscEqflyggcegnannfytweacddacwri		SEQ ID NO. 71
Human TFPI-2 DOMAIN 2	VPKVCRLQV- SVDDQCEGSTEKYFFNLSMTCEKFFSGGCHRNR- IENRFPDEATCMGFCAPK		SEQ ID NO. 72
DPI-3.2.1	vpkvcrlqvETGpCRgKfkyffnlsmtceTfvYggcEGnrnrfpdeatcmgfcapk		SEQ ID NO. 73
Human TFPI-2 DOMAIN 3	IPSPCYSPKDEGLCSANVTRYFYNPRYRTCDAFYTGCGGNDNNFVSREDCKRACAKA		SEQ ID NO. 74
DPI-3.3.1	ipsfcyspkdTgPcraRftryfynpryrtcdafyGgcegnndnnfvsredckracaka		SEQ ID NO. 75
HUMAN ITI-K1	KEDSCQLGYSAGPCMGMTSRIFYNGTSMACETFQYGGCMGNNFVTEKECLQTCRTV		SEQ ID NO. 76
DPI-4.1.1	kedscqlgyEagpcRgKfryfyngtSmacetfVyggeGgnnnfvtekeclqtcrtv		SEQ ID NO. 77
Human ITI-K2	TVAACNLPiVRGPCRAFQLWAFDAVKGKCVLFPYGGCQGNKNFYSEKECREYCGVP		SEQ ID NO. 78
DPI-4.2.1	tvaacnlpIDTgpcraRfqlwafdavkgkcvlFvyggcqqgnknfysekecreycgvp		SEQ ID NO. 79
PROTEASE NEXIN-II	VREVCSEQAETGPCRAMISRWFYDVTGKCAPFFYGGCGGNNRNFDTTEEYCMAVCGSA		SEQ ID NO. 80
DPI-5.1	vrevcseqaetgpcraRfswfyfdvtgkcapffYggcgggnrnnfdteeycmavcgsa		SEQ ID NO. 81
DPI-5.2	vrevcseqaetgpcraRfswfyfdvtgkcePfiyggcgggnrnnfdteeycmavcgsa		SEQ ID NO. 82
HKI B9	LPNVCAFPMKGCQTYMTRWFFNFETGEBELFAYGGCGGNSNNFLRKEKCEKFCCKFT		SEQ ID NO. 124
DPI-6.1	lpnvcafpmETgpcRARftrwffnfetgecelfayggcggnsnnflrkekeckfckft		SEQ ID NO. 125
DPI-6.2	lpnvcafpmETgpcRARFDrffnfetgecelfVyggeggnsnnflrkekeckfckft		SEQ ID NO. 126
DPI-4.2.2	tvaacnlpivTgpcraRfqlwafdavkgkcvlFpyggcqqgnknfysekecreycgvp		SEQ ID NO. 130

TABLE 4-continued

<u>Kunitz domains, some of which inhibit plasmin</u>			
Protein identifier	Amino-acid Sequence	Affinity K_D	SEQ ID NO.
DPI-4.2.3	tvaacnlpivTgpcraRFqRwafdavkgkcvlfpyggcqqgngnkfysekecreycgvp		SEQ ID NO. 131
DPI-4.2.4	tvaacnlpivTgpcraRFqRwafdavkgkcvlfVyggcqqgngnkfysekecreycgvp		SEQ ID NO. 132
DPI-4.2.5	tvaacnlpieTgpcraRfDRwafdavkgkceTfVyggcGgngnkfysekecreycgvp		SEQ ID NO. 133
DPI-7.1	rpdfcleppytgpcRarFiryfynakaglcqtfvyygcrakrnnfksaedcmrtcgga		SEQ ID NO. 134
DPI-7.2	rpdfcleppytgpcRarFiryfynakaglcqtfvyygcrakrnnfksaedcmrtcgga		SEQ ID NO. 135
DPI-7.3	rpdfcleppytgpcRarFiryfynakaglcqtfvyygGakrnnfksaedcmrtcgga		SEQ ID NO. 136
DPI-7.4	rpdfcleppDtgpCRarFDryfynakaglcEtfvyygGakrnnfksaedcmrtcgga	⊙	SEQ ID NO. 137

Under "Affinity", "(A)" means the K_D is likely to be less than that of BPTI (viz.300 pM), "(B)" means K_D is likely to be less than 2 nM, and "(C)" means that K_D is likely to be less than 20 nM.

TABLE 5

vgDNA for LACI-DI to vary residues 10, 11, 13, 15, 16, 17, & 19 for plasmin in view of App-I (now known not to be very potent)

										N K
M	H	S	F	C	A	F	K	A	D E	
1	2	3	4	5	6	7	8	9	10	
5'- cctcct	atgcat	tcc	ttc	tgc	gcc	ttc	aag	gct	RaS	-
	<u>NsiI</u>									
										C W
										F Y
										L P
										Q H
										M I
										N T
N S					T					S K
I T					N K					V R
A G	A P				I M					D A
D V	G S	T C	K R	A G	R S	F I	E G	R		
11	12	13	14	15	16	17	18	19	20	
Rnt	ggt	Nct	tgt	aRa	gSt	aNS	wtc	NNS	cgt	
										F C
L W	F	F	N	I	F	T	R			
21	22	23	24	25	26	27	28			
tKS	ttc	ttc	aac	atc	ttc	acg	cgT	tcctcc-3'	(SEQ ID NO. 83)	
3'-g	aag	ttg	tag	aag	tgc	gca	agggagg-5'	(SEQ ID NO. 84)		
Tm - 80° C.					<u>MluI</u>					

DNA: 262,144 * 4 = 1,048,576
 protein: 143,360 * 4 = 573,440
 The amino acid seq has SEQ ID NO. 85.
 This variegation allows the AppI sequence to appear in the P6-P6' positions.

TABLE 6

<u>LACI-K1 derivatives selected for Plasmin binding</u>				
Ident	111111111122	Phage	K_D	SEQ ID NO.
	012345678901	DIFFS	Binding	(pM)
Consensus	ETGPCRARFERW	0		SEQ ID NO. 88
LACI-K1	ddgpckaimkrf	7		SEQ ID NO. 2

TABLE 6-continued

<u>LACI-K1 derivatives selected for Plasmin binding</u>					
Ident	111111111122 012345678901	DIFFS	Phage Binding	K_D (pM)	SEQ ID NO.
SPI31	-----G--	1			SEQ ID NO. 89
SPI11	-----D--	1	3.2 x	88	SEQ ID NO. 40
SPI15	-S-----D--	2	2.5 x		SEQ ID NO. 90
SPI24	D----G----L	3			SEQ ID NO. 91
SPI33	---S--G--D--	3			SEQ ID NO. 92
SPI34	-V-----S-P--	3			SEQ ID NO. 93
SPI26	-----T-P-F	3			SEQ ID NO. 94
SPI37	-V-----S-H--	3			SEQ ID NO. 95
SPI32	D-----S-G--	3			SEQ ID NO. 96
SPI12	-----GM-P--	3			SEQ ID NO. 97
SPI36	-G-----N-F	3			SEQ ID NO. 98
SPI08	DG-----F	3	2.6 x		SEQ ID NO. 42
SPI38	-----IS-F	3			SEQ ID NO. 99
SPI18	-G-----K---F	3			SEQ ID NO. 100
SPI23	-G-----K-Q--	3	1.25 x		SEQ ID NO. 43
SPI35	DS-A--G-----	4			SEQ ID NO. 101
SPI02	DS---G---F	4	0.83 x		SEQ ID NO. 102
SPI25	D-----S-P-L	4			SEQ ID NO. 103
SPI17	-V-----IQ-F	4	0.09 x		SEQ ID NO. 104
SPI05	-S-----K-A-F	4	0.64 x		SEQ ID NO. 105
SPI13	-G-----K-A-F	4			SEQ ID NO. 106
SPI07	D--S---KI---	4			SEQ ID NO. 107
SPI03	DS---K---D--	4	0.48 x		SEQ ID NO. 108
SPI06	DG---KG-----	4			SEQ ID NO. 109
SPI16	-V-A-KG--H--	5	0.22 x		SEQ ID NO. 110
SPI04	DG-----S-P-F	5			SEQ ID NO. 111
SPI01	DS-A---M-H-F	6	0.25 x		SEQ ID NO. 112
SPI14	DS-A---K-R--	5			SEQ ID NO. 113
SPI28	DS-T-K---P-F	6			SEQ ID NO. 114
SPI27	-----KGKIA-F	6			SEQ ID NO. 115
SPI21	DS-A-KGK----	6	0.38 x		SEQ ID NO. 116
SPI22	DG---KGK-P-F	7	2.0 x		SEQ ID NO. 44

"Diffs" is the number of differences from the Consensus.
 "Phage Binding" is the binding of phage that display the
 named protein relative to binding of phage that display BPTI.

TABLE 7

Variation of Residues 31, 32, 34, and 39													
5'-cctcct acg cgt cag tgc -													
<u>MluI</u>													
F S F S													
Y C Y C													
L P L P													
H R H R	G												
W I W I	N D												
T M T M	H R												
N K N K	Y C												
V A V A	A V												
D G D G	I T												
E Q E Q	S P C												
31	32	33	34	35	36	37	38	39	40	41	42		
NNS NNS ttc NNT tRS ggt ggt tgt RRg ggt aac cag -													
<u>BstEII</u>													

gtcgtgctcttttagcacgacctg-3' (SEQ ID NO. 86)

The amino acid sequence has SEQ ID NO. 87.
 The EcoRI site is erased; thus, cleavage with EcoRI can be used to eliminate (or at least greatly reduce) parental DNA.
 There are 262,144 DNA sequences and 72,000 protein sequences.

TABLE 8

Selectants for plasmin binding with variegation of second loop						
111111111122		3333333334		# Diffs		
Id	012345678901	1234567890	C1	C	K1	
Con1	ETgPcRAKFDrW	EafVYggcGg			10	SEQ ID NO. 45
SPI47	DA-----E-F	-E-I----E-	7	(5)	5	SEQ ID NO. 57
SPI51	DA-----E-F	T--F-----	6	(4)	9	SEQ ID NO. 54
SPI52	DA-----E-F	---L-----	5	(3)	8	SEQ ID NO. 52
SPI46	DV-----E-F	---L----E-	6	(3)	7	SEQ ID NO. 53
SPI41	DA-----R-E-F	DT-L----E-	9	(6)	8	SEQ ID NO. 59
SPI42	-----G-----	QG-----	3	(3)	12	SEQ ID NO. 48
SPI43	-----G-----	-S-H----K-	4	(4)	11	SEQ ID NO. 51
SPI56	-----G-----	AQ-----E-	4	(3)	11	SEQ ID NO. 50
SPI59	-----	NT-----	2	(2)	11	SEQ ID NO. 47
SPI60	-----	-P-----E-	2	(1)	9	SEQ ID NO. 46
SPI55	-V-----	HL-T-----	4	(4)	11	SEQ ID NO. 49
SPI49	-----IP-L	-P-IW-----	6	(6)	10	SEQ ID NO. 56
SPI57	DS---KGR-G-L	T--DW-----	10	(8)	11	SEQ ID NO. 60
SPI53	-----KGS----	NV-R----R-	7	(7)	11	SEQ ID NO. 58
SPI54	DS-----R----	TR-P-----	6	(4)	10	SEQ ID NO. 55

TABLE 8-continued

Selectants for plasmin binding with variegation of second loop					
	111111111122	3333333334	# Diffs		
Id	012345678901	1234567890	C1	C	K1
SPI11	-----R----	eefiyggceg	(1)	(4)	7 SEQ ID NO. 40
LAC11	ddgpckaimkrf	eefiyggceg	10	(7)	0 SEQ ID NO. 2

See notes below.

In the Table, "-" means that the protein has the consensus (Con1) type. Con1 contains the most common type at each position; amino acids shown in Con1 were not varied. Four positions (10, 31, 34, and 39) showed significant toleration for a second type, leading to 15 subsidiary consensus sequences: Con2-Con16. The column "# Diffs" shows the number of differences from CON1 under "C1", the differences with the closest of Con1-Con16 under "C", and the differences from LAC1-K1 under "K1". SPI11 was selected from a library in which residues 31-39 were locked at the wild-type.

SPI11 <BPTI <SPI23 ≈ SPI51 <SPI47 <QS4 <SPI22 <SPI54 <SPI43
Highly very potent
Superior potent

TABLE 9

Conservative and Semiconservative substitutions			
Initial AA type	Category	Conservative substitution	Semi-conservative substitution
A	Small non-polar or slightly polar	G, S, T	N, V, P, (C)
C	free SH disulfide	A, M, L, V, I	F, G
D	acidic, hydrophilic	E, N, S, T, Q	K, R, H, A
E	acidic, hydrophilic	D, Q, S, T, N	K, R, H, A
F	aromatic	W, Y, H, L, M	I, V, (C)
G	Gly-only conformation "normal"	nothing	nothing
	conformation	A, S, N, T	D, E, H, I, K, L, M, Q, R, V
H	amphoteric aromatic	Y, F, K, °R	L, M, A, (C)
I	aliphatic, branched β carbon	V, L, M, A	F, Y, W, G (C)
K	basic	R, H	Q, N, S, T, D, E, A
L	aliphatic	M, I, V, A	F, Y, W, H, (C)
M	hydrophobic	L, I, V, A	Q, F, Y, W, (C), (R), (K), (E)
N	non-polar hydrophilic	S, T, (D), Q, A, G, (E)	K, R

TABLE 9-continued

Conservative and Semiconservative substitutions			
Initial AA type	Category	Conservative substitution	Semi-conservative substitution
P	inflexible	V, I	A, (C), (D), (E), F, H, (K), L, H, N, Q, (R), S, T, W, Y
Q	aliphatic plus amide	N, E, A, S, T, D	M, L, K, R
R	basic	K, Q, H	S, T, E, D, A,
S	hydrophilic	A, T, G, N	D, E, R, K
T	hydrophilic	A, S, G, N, V	D, E, R, K, I
V	aliphatic, branched β carbon	I, L, M, A, T	P, (C)
W	aromatic	F, Y, H	L, M, I, V, (C)
Y	aromatic	F, W, H	L, M, I, V, (C)

[0152] Changing from A, F, H, I, L, M, P, V, W, or Y to C is semiconservative if the new cysteine remains as a free thiol.

[0153] Changing from M to E, R, K is semiconservative if the ionic tip of the new side group can reach the protein surface while the methylene groups make hydrophobic contacts.

[0154] Changing from P to one of K, R, E, or D is semiconservative if the side group is on or near the surface of the protein.

TABLE 10

Plasmin-inhibiting Kunitz domain derivatives of LAC1-K1									
		Consensus #1		Consensus #2		Consensus #3		Consensus #4	
Position	Type	Status	Type	Status	Type	Status	Type	Status	
10	D	fixed	D	fixed	E/D	S-S	D/E	S-S	
11	D	fixed	D	fixed	T/S	G-S	T/A	G-S	
12	G	fixed	G	fixed	G	fixed	G	fixed	
13	P	Abs-S	P	VS-S	P	VS-S	P	Abs-S	

TABLE 10-continued

<u>Plasmin-inhibiting Kunitz domain derivatives of LACI-K1</u>								
	<u>Consensus #1</u>		<u>Consensus #2</u>		<u>Consensus #3</u>		<u>Consensus #4</u>	
Position	Type	Status	Type	Status	Type	Status	Type	Status
14	C	fixed	C	fixed	C	fixed	C	fixed
15	K	fixed	K	fixed	R	S-S	R	S-S
16	A	Abs-S	A	Abs-S	A	VS-S	A	S-S
17	R	Abs-S	R	VS-S	R/K	S-S	K	S-S
18	F	Abs-S	F	Abs-S	F	VS-S	F	VS-S
19	E	Abs-S	E	Abs-S	E/P/D	S-S	D/E	VS-S
20	R	fixed	R	fixed	R	fixed	R	fixed
21	F	fixed	F	fixed	W/F	weak-Sel	W/F	weak-Sel
31	E	S-S	E	S-S	E	fixed	E/t	G-S
32	Q	G-S	Q	G-S	E	fixed	A/T	Strong for no charge, weak for type fixed
33	F	fixed	F	fixed	F	fixed	F	fixed
34	—	no consensus	T/S	weak	I	fixed	V/L/I	Weak
35	Y	fixed	Y	fixed	Y	fixed	Y	S-S
39	—	no consensus	G	weak Sel.	E	fixed	G/E	some-Sel.

Abs-S Absolute Selection
 VS-S Very Strong Selection
 S-S Strong Selection
 G-S Good Selection

TABLE 11

<u>High Specificity Designed Plasmin Inhibitors</u>		
Ident	Sequence	SEQ ID NO.
SPI11	mhsfcafkaETgPcRARFDrWffniftrqceefiyggcegnqnrfsleackmctrd	SEQ ID NO. 40
SPI11-R15A	mhsfcafkaETgPcAARFDrWffniftrqceefiyggcegnqnrfsleackmctrd	SEQ ID NO. 117
SPI1-R15G	mhsfcafkaETgPcGARFDrWffniftrqceefiyggcegnqnrfsleackmctrd	SEQ ID NO. 118
SPI11-R15N-E32A	mhsfcafkaETgPcNARFDrWffniftrqceAfiyggcegnqnrfsleackmctrd	SEQ ID NO. 117

TABLE 12

vgDNA for LACI-D1 to vary residues 10, 11, 12, 13, 14, 15, 16, 17, 19, 20, 21, 37, 38, and 39 for plasmin in view of App-I and SPI11

	M	H	S	F	C	A	F	K	A	D E
5'-	1	2	3	4	5	6	7	8	9	10
cctcct	atgcat	tcc	ttc	tgc	gcc	ttc	aag	gct	GaS	
	NsiI									
	Y	L	F	S	F	S				
	F	S	Y	C	Y	C				
	C	P	L	P	L	P				
	H	R	H	R	H	R				
	I	T	I	T	I	T	N S			
	N	V	N	V	N	V	I	T	V	D
S P	A	D	A	D	A	D	V	D	A	E
T	A	G	G	G	K	R	A	G	R	Q
11	12	13	14	15	16	17	18	19	20	
NCT	NNT	NNT	NNT	aRa	RNT	aRa	ttc	gNS	cRt	

TABLE 12-continued

vgDNA for LACI-D1 to vary residues 10, 11, 12, 13,
14, 15, 16, 17, 19, 20, 21, 37, 38, and 39 for
plasmin in view of App-I and SPII1

```

F|C
L|W F F N I F T R Q C
21 22 23 24 25 26 27 28 29 30
|tKS|ttc|ttc|aac|atc|ttc|acg cgt|cag|tgc|-3'
3'-|aag aag ttg tag aag tgc gca gtc acg-
      |MluI|

      F|S F|S
      Y|C Y|L
      L|P P|H I|M
      H|R R|I T|N
      I|T T|N K|S
      N|V V|A R|V
      A|D D|G A|E

E A F V Y G G C G D G N Q
31 32 33 34 35 36 37 38 39 40 41 42
|cag|cct|ttc|ctt|tac|ggc|nnc|nnc|rns|ggc|aac|cag|
ctc cga aag caa atg cca nna nna yns cca ttg gtc cctcctcc-5'
      |BstEII|

```

First (top) strand of DNA has SEQ ID NO. 120.

Second (bottom) strand of DNA has SEQ ID NO. 121.

The amino-acid sequence has SEQ ID NO. 122.

The top strand for codons 31-42 (shown stricken) need not be synthesized, but is produced by PCR from the strands shown.

There are 1.37×10^{11} DNA sequences that encode 4.66×10^{10} amino-acid sequences.

TABLE 14

Definition of a Kunitz Domain (SEQ ID NO. 123)

```

1      2      3      4      5
1234567890123456789012345678901234567890123456789012345678
xxxxCxXXXXXXXXGxCXXXXXXXXXXXXXXXxCxxFxXXGCxXxxXxXXXXXXXXCxXxx

```

Where:

X1, X2, X3, X4, X58, X57, and X56 may be absent,

X21 = Phe, Tyr, Trp,

X22 = Tyr or Phe,

X23 = Tyr or Phe,

X35 = Tyr or Trp,

X36 = Gly or Ser,

X40 = Gly or Ala,

X43 = Asn or Gly, and

X45 = Phe or Tyr

TABLE 15

Substitutions to confer high affinity for plasmin on KuDoms

Position	Allowed types
10	Asp, Glu, Tyr
11	Thr, Ala, Ser, Val, Asp
12	Gly
13	Pro, Leu, Ala
14	Cys

TABLE 15-continued

Substitutions to confer high affinity for plasmin on KuDoms

Position	Allowed types
15	Arg, Lys
16	Ala, Gly
17	Arg, Lys, Ser
18	Phe, Ile
19	Glu, Asp, Pro, Gly, Ser, Ile

TABLE 15-continued

<u>Substitutions to confer high affinity for plasmin on KuDoms</u>	
Position	Allowed types
20	Arg
21	Phe, Trp, Tyr
31	Asp, Glu, Thr, Val, Gln, Ala
32	Thr, Ala, Glu, Pro, Gln
34	Val, Ile, Thr, Leu, Phe, Tyr, His, Asp, Ala, Ser
35	Tyr, Trp
36	Gly
37	Gly
38	Cys
39	Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys, Met

TABLE 16

<u>Summary of Sequences selected from First LACI-K1 library for binding to Plasmin</u>			
BPTI # (BPTI type)	(LACI-K1)	Residues Allowed in Library	Preferred Residues
13 (P)	P	LHPR	PL
16 (A)	A	AG	AG
17 (R)	I	FYLHINA SCPRTVDG	RS
18 (I)	M	all	F
19 (I)	K	LWQMKAG SPRTVE	EQ
31 (Q)	E	EQ	EQ
32 (T)	E	EQ	QE
34 (V)	I	all	TYHDRVILSF
39 (R)	E	all	GADRQFEMLVKNH

In Table 15 the bold residue types are preferred

TABLE 17

<u>Distribution of sequences selected from first library:</u>																				
Position	A	C	D	E	F	G	H	I	K	L	M	N	P	Q	R	S	T	V	W	Y
13	x	x	x	x	x	x	0	x	x	1	x	x	31*	x	0	x	x	x	x	x
16	31*	x	x	x	x	1	x	x	x	x	x	x	x	x	x	x	x	x	x	x
17	0	0	0	x	0	0	0	0*	x	0	x	0	0	x	30	2	0	0	x	0
18	0	0	0	0	32	0	0	0	0	0	0*	0	0	0	0	0	0	0	0	0
19	0	x	x	31	x	0	x	x	0*	0	0	x	0	1	0	0	0	0	0	x
31	x	x	x	28*	x	x	x	x	x	x	x	x	x	4	x	x	x	x	x	x
32	x	x	x	9*	x	x	x	x	x	x	x	x	x	23	x	x	x	x	x	x
34	2	0	1	0	1	0	2	5*	0	2	0	0	0	0	1	3	8	5	0	2
39	3	0	3	2*	1	10	1	0	2	2	2	1	0	3	1	0	0	1	0	0

TABLE 18

<u>Distribution of amino-acid types at varied residues in proteins selected for plasmin binding from third library.</u>																				
Position	A	C	D	E	F	G	H	I	K	L	M	N	P	Q	R	S	T	V	W	Y
10	x	x	7*	8	x	x	x	x	0	x	x	0	x	x	x	x	x	x	x	x
11	4	x	0*	x	x	0	x	0	x	x	x	0	x	x	x	2	7	2	x	x
13	0	x	x	x	x	x	x	x	x	x	x	x	15*	x	x	0	0	x	x	x
15	x	x	x	x	x	x	x	x	2*	x	x	x	x	x	13	x	x	x	x	x
16	10*	x	x	x	x	5	x	x	x	x	x	x	x	x	x	x	x	x	x	x
17	x	x	x	x	x	x	x	0*	11	x	0	0	x	x	3	1	0	x	x	x
18	x	x	x	x	14	x	x	1	x	x	x*	x	x	x	x	x	x	x	x	x
19	0	0	8	5	0	1	0	0	0*	0	0	0	1	0	0	0	0	0	0	0
21	x	0	x	x	5*	x	x	x	x	2	x	x	x	x	x	x	x	x	8	x
31	1	0	1	6*	0	0	1	0	0	0	0	2	0	1	0	0	3	0	0	0
32	4	0	0	1*	0	1	0	0	0	1	0	0	2	1	1	1	2	1	0	0
34	0	0	1	x	1	0	1	2*	x	3	x	0	1	x	1	0	1	4	x	0
35	x	0	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	2	13*
39	x	x	x	5*	x	8	x	x	1	x	x	x	x	x	1	x	x	x	x	x

TABLE 23

Specificity Results Obtained with KuDoms Displayed on gIIIp of M13					
KuDom Displayed	Target				
	Plasmin	Thrombin	Kallikrein	Trypsin	Trypsin, 2 washes
LACI-K1	1.0	1.0	1.0	1.0	1.0
QS4	52.	0.7	0.9	4.5	0.5
BPTI	88.	1.1	1.7	0.3	0.8

LACI-K1 phage for each Target was taken as unit binding and the other display phage are shown as relative binding. BPTI::III phage are not easily liberated from trypsin.

CITATIONS

- [0155] ADEL86: Adelman et al., *Blood* (1986) 68(6)1280-1284.
 [0156] ANBA88: Anba et al., *Biochimie* (1988) 70(6)727-733.
 [0157] AUER88: Auerswald et al., *Bio Chem Hoppe-Seyler* (1988), 369 (Supplement):27-35.
 [0158] BANE90: Baneyx & Georgiou, *J Bacteriol* (1990) 172(1)491-494.
 [0159] BANE91: Baneyx & Georgiou, *J Bacteriol* (1991) 173(8)2696-2703.
 [0160] BROW91: Browne et al., *GeneBank* entry M74220.
 [0161] BROZ90: Broze et al., *Biochemistry* (1990) 29:7539-7546.

TABLE 24

Mat α <i>S. cerevisiae</i> expression vectors:	
Mat α 1 (Mfa8)	
5' - ... AAA AGG CCT CGA G... - 3'	
StuI XhoI	
Mat α 2 (after introduction of a linker into StuI-cut DNA)	
5' AAA AGG GAA GCG GCC GAG CCA TGG GGC GCC TAA TAG CTC GAG 3'	
EaqI StuI Kasi XhoI	
Mat α -LACI-K1	
5' AAA AGG GAA GCG GCC GAG atg cat tcc ttc tgc gct ttc aaa	
EaqI NsiI	
9 10 11 12 13 14 15 16 17 18 19 20	
A D D G P C K A I M K R	
gct gat gac ggt ccg tgt aaa gct atc atg aaa cgt	
RsrII BspHI	
21 22 23 24 25 26 27 28 29 30	
F F F N I F T R Q C	
ttc ttc ttc aac att ttc acg cgt cag tgc	
MluI	
31 32 33 34 35 36 37 38 39 40 41 42	
E E F I Y G G C E G N Q	
gag gaA ttc att tac ggt ggt tgt gaa ggt aac cag	
EcoRI BstEII	
43 44 45 46 47 48 49 50	
N R F E S L E E	
aac cgG ttc gaa tct ctA gag gaa	
AqeI BstBI XbaI	
51 52 53 54 55 56 57 58 59 60	
C K K M C T R D G A	
tgt aag aag atg tgc act cgt gac ggc gcc TAA TAG CTC GAG - 3'	
Kasi XhoI	

We expect that Mat α pre sequence is cleaved before GLU_a -ALA_a-

- [0162] COLM87: Colman et al., Editors, Hemostasis and Thrombosis, Second Edition, 1987, J. B. Lippincott Company, Philadelphia, Pa.
- [0163] DENN94a: Dennis & Lazarus, *J Biological Chem* (1994) 269:22129-22136.
- [0164] DENN94b: Dennis & Lazarus, *J Biological Chem* (1994) 269:22137-22144.
- [0165] EIGE90: Eigenbrot et al., *Protein Engineering* (1990), 3(7):591-598.
- [0166] ELLI92: Ellis et al., *Ann NY Acad Sci* (1992) 667: 13-31.
- [0167] FIDL94: Fidler & Ellis, *Cell* (1994) 79:185-188.
- [0168] FRAE89: Fraedrich et al., *Thorac Cardiovasc Surg* (1989) 37(2):89-91.
- [0169] GARD93: Gardell, *Toxicol Pathol* (1993) 21(2):190-8.
- [0170] GIRA89: Girard et al., *Nature* (1989), 338:518-20.
- [0171] GIRA91: Girard et al., *J. BIOL. CHEM.* (1991) 266: 5036-5041.
- [0172] GREG93: Gregg et al., *Bio/Technology* (1993)11: 905-910.
- [0173] HOOV93: Hoover et al., *Biochemistry* (1993) 32:10936-43.
- [0174] HYNE90: Hynes et al., *Biochemistry* (1990), 29:10018-10022.
- [0175] KIDO88: Kido et al., *J Biol Chem* (1988), 263: 18104-7.
- [0176] KIDO90: Kido et al., *Biochem & Biophys Res Comm* (1990), 167(2):716-21.
- [0177] KURJ82: Kurjan and Herskowitz, *Cell* (1982) 30:933-943.
- [0178] LASK80: Laskowski & Kato, *Ann Rev Biochem* (1980), 49:593-626.
- [0179] LEAT91: Leatherbarrow & Salacinski, *Biochemistry* (1991) 30(44):10717-21.
- [0180] LOHM93: Lohmann & J Marshall, *Refract Corneal Surg* (1993) 9(4):300-2.
- [0181] LUCA83: Lucas et al., *J Biological Chem* (1983) 258(7):4249-56.
- [0182] MANN87: Mann & Foster, Chapter 10 of COLM87.
- [0183] MIYA85: Miyajima et al., *Gene* (1985) 37:155-161.
- [0184] NEUH89: Neuhaus et al., *Lancet* (1989) 2(8668) 924-5.
- [0185] NOVO89: Novotny et al. *J. BIOL. CHEM.* (1989) 264:18832-18837.
- [0186] OREI94: O'Reilly et al., *Cell* (1994) 79:315-328.
- [0187] PARK86: Park & Tulinsky, *Biochemistry* (1986) 25(14):3977-3982.
- [0188] PUTT89: Putterman, *Acta Chir Scand* (1989) 155 (6-7):367.
- [0189] ROBB87: Robbins, Chapter 21 of COLM87
- [0190] SCHE67: Schechter & Berger, *Biochem Biophys Res Commun* (1967) 27:157-162.
- [0191] SCHE68: Schechter & Berger, *Biochem Biophys Res Commun* (1969) 32:898-902.
- [0192] SCHN86: Schnabel et al., *Biol Chem Hoppe-Seyler* (1986), 367:1167-76.
- [0193] SHER89: Sheridan et al., *Dis Colon Rectum* (1989) 32(6):505-8.
- [0194] SPRE94: Sprecher et al., *Proc Natl Acad Sci USA* 91:3353-3357 (1994)
- [0195] VAND92: van Dijn et al., *EMBO J* (1992) 11(8) 2819-2828.
- [0196] VARA83: Varadi & Patthy, *Biochemistry* (1983) 22:2440-2446.
- [0197] VARA84: Varadi & Patthy, *Biochemistry* (1984) 23:2108-2112.
- [0198] VEDV91: Vedvick et al., *J Ind Microbiol* (1991) 7:197-201.
- [0199] WAGN92: Wagner et al., *Biochem Biophys Res Comm* (1992) 186:1138-1145.
- [0200] WUNT88: Wun et al., *J. BIOL. CHEM.* (1988) 263: 6001-6004.

SEQUENCE LISTING

<160> NUMBER OF SEQ ID NOS: 153

<210> SEQ ID NO 1

<211> LENGTH: 304

<212> TYPE: PRT

<213> ORGANISM: Homo sapiens

<400> SEQUENCE: 1

Met Ile Tyr Thr Met Lys Lys Val His Ala Leu Trp Ala Ser Val Cys
1 5 10 15

Leu Leu Leu Asn Leu Ala Pro Ala Pro Leu Asn Ala Asp Ser Glu Glu
20 25 30

Asp Glu Glu His Thr Ile Ile Thr Asp Thr Glu Leu Pro Pro Leu Lys
35 40 45

Leu Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys
50 55 60

Ala Ile Met Lys Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu
65 70 75 80

Glu Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser
85 90 95

-continued

Leu Glu Glu Cys Lys Lys Met Cys Thr Arg Asp Asn Ala Asn Arg Ile
 100 105 110
 Ile Lys Thr Thr Leu Gln Gln Glu Lys Pro Asp Phe Cys Phe Leu Glu
 115 120 125
 Glu Asp Pro Gly Ile Cys Arg Gly Tyr Ile Thr Arg Tyr Phe Tyr Asn
 130 135 140
 Asn Gln Thr Lys Gln Cys Glu Arg Phe Lys Tyr Gly Gly Cys Leu Gly
 145 150 155 160
 Asn Met Asn Asn Phe Glu Thr Leu Glu Glu Cys Lys Asn Ile Cys Glu
 165 170 175
 Asp Gly Pro Asn Gly Phe Gln Val Asp Asn Tyr Gly Thr Gln Leu Asn
 180 185 190
 Ala Val Asn Asn Ser Leu Thr Pro Gln Ser Thr Lys Val Pro Ser Leu
 195 200 205
 Phe Glu Phe His Gly Pro Ser Trp Cys Leu Thr Pro Ala Asp Arg Gly
 210 215 220
 Leu Cys Arg Ala Asn Glu Asn Arg Phe Tyr Tyr Asn Ser Val Ile Gly
 225 230 235 240
 Lys Cys Arg Pro Phe Lys Tyr Ser Gly Cys Gly Gly Asn Glu Asn Asn
 245 250 255
 Phe Thr Ser Lys Gln Glu Cys Leu Arg Ala Cys Lys Lys Gly Phe Ile
 260 265 270
 Gln Arg Ile Ser Lys Gly Gly Leu Ile Lys Thr Lys Arg Lys Arg Lys
 275 280 285
 Lys Gln Arg Val Lys Ile Ala Tyr Glu Glu Ile Phe Val Lys Asn Met
 290 295 300

<210> SEQ ID NO 2
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Homo sapiens

<400> SEQUENCE: 2

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15
 Ile Met Lys Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30
 Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45
 Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 3
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 3

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15
 Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

-continued

Phe Thr Tyr Gly Gly Cys Arg Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 4
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 4

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Tyr Tyr Gly Gly Cys Asp Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 5
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 5

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe His Tyr Gly Gly Cys Asp Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 6
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 6

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Asp Tyr Gly Gly Cys Ala Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 7
 <211> LENGTH: 58
 <212> TYPE: PRT

-continued

```

<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 7
Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1             5             10             15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Gln Glu
 20             25             30
Phe Arg Tyr Gly Gly Cys Asp Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 8
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 8
Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1             5             10             15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Gln Gln
 20             25             30
Phe Tyr Tyr Gly Gly Cys Gln Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 9
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 9
Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1             5             10             15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20             25             30
Phe Ala Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 10
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 10
Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1             5             10             15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Gln Gln

```

-continued

20	25	30
Phe Val Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu		
35	40	45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp		
50	55	

<210> SEQ ID NO 11
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 11

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala		
1	5	10
15		
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln		
20	25	30
Phe Thr Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu		
35	40	45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp		
50	55	

<210> SEQ ID NO 12
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 12

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala		
1	5	10
15		
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu		
20	25	30
Phe Thr Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu		
35	40	45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp		
50	55	

<210> SEQ ID NO 13
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 13

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala		
1	5	10
15		
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln		
20	25	30
Phe Ile Tyr Gly Gly Cys Gln Gly Asn Gln Asn Arg Phe Glu Ser Leu		
35	40	45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp		
50	55	

<210> SEQ ID NO 14

-continued

<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 14

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
20 25 30

Phe Ile Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 15
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 15

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
20 25 30

Phe Ile Tyr Gly Gly Cys Phe Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 16
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 16

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
20 25 30

Phe His Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 17
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 17

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
1 5 10 15

-continued

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Val Tyr Gly Gly Cys Ala Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 18
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 18

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Leu Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 19
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 19

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Ile Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 20
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 20

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Val Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

-continued

<210> SEQ ID NO 21
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 21

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
1 5 10 15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Val Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 22
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 22

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Leu Cys Lys Gly
1 5 10 15
Arg Phe Gln Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 23
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 23

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
1 5 10 15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
20 25 30
Phe Thr Tyr Gly Gly Cys Met Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 24
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 24

-continued

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1           5           10           15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
      20           25           30
Phe Ser Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
      50           55

```

```

<210> SEQ ID NO 25
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 25

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1           5           10           15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
      20           25           30
Phe Leu Tyr Gly Gly Cys Leu Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
      50           55

```

```

<210> SEQ ID NO 26
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 26

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1           5           10           15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
      20           25           30
Phe Ser Tyr Gly Gly Cys Gln Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
      50           55

```

```

<210> SEQ ID NO 27
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 27

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1           5           10           15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
      20           25           30
Phe Ala Tyr Gly Gly Cys Ala Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp

```


-continued

50 55

<210> SEQ ID NO 28
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 28

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Ile Tyr Gly Gly Cys Val Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 29
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 29

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ser Tyr Gly Gly Cys Lys Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 30
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 30

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Val Tyr Gly Gly Cys Lys Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 31
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

-continued

<400> SEQUENCE: 31

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1             5             10             15
Ser Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
      20             25             30
Phe Thr Tyr Gly Gly Cys Asn Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

<210> SEQ ID NO 32

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 32

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1             5             10             15
Ser Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
      20             25             30
Phe Thr Tyr Gly Gly Cys Leu Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

<210> SEQ ID NO 33

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 33

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1             5             10             15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
      20             25             30
Phe Phe Tyr Gly Gly Cys His Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

<210> SEQ ID NO 34

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 34

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1             5             10             15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
      20             25             30
Phe Thr Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35             40             45

```

-continued

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 35
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 35

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Thr Tyr Gly Gly Cys Met Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 36
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 36

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Thr Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 37
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 37

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Gln
 20 25 30

Phe Val Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 38
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Bos taurus

-continued

<400> SEQUENCE: 38

Arg Pro Asp Phe Cys Leu Glu Pro Pro Tyr Thr Gly Pro Cys Lys Ala
 1 5 10 15

Arg Ile Ile Arg Tyr Phe Tyr Asn Ala Lys Ala Gly Leu Cys Gln Thr
 20 25 30

Phe Val Tyr Gly Gly Cys Arg Ala Lys Arg Asn Asn Phe Lys Ser Ala
 35 40 45

Glu Asp Cys Met Arg Thr Cys Gly Gly Ala
 50 55

<210> SEQ ID NO 39

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Homo sapiens

<400> SEQUENCE: 39

Val Arg Glu Val Cys Ser Glu Gln Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Met Ile Ser Arg Trp Tyr Phe Asp Val Thr Glu Gly Lys Cys Ala Pro
 20 25 30

Phe Phe Tyr Gly Gly Cys Gly Glu Asn Arg Asn Asn Phe Asp Thr Glu
 35 40 45

Glu Tyr Cys Met Ala Val Cys Gly Ser Ala
 50 55

<210> SEQ ID NO 40

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 40

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 41

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 41

Met His Ser Phe Cys Ala Phe Lys Ala Glu Ser Gly Pro Cys Arg Ala
 1 5 10 15

Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp

-continued

```

50                55

<210> SEQ ID NO 42
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 42
Met His Ser Phe Cys Ala Phe Lys Ala Asp Gly Gly Pro Cys Arg Ala
 1                5                10                15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
                20                25                30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
                35                40                45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50                55

<210> SEQ ID NO 43
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 43
Met His Ser Phe Cys Ala Phe Lys Ala Glu Gly Gly Pro Cys Arg Ala
 1                5                10                15
Lys Phe Gln Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
                20                25                30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
                35                40                45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50                55

<210> SEQ ID NO 44
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 44
Met His Ser Phe Cys Ala Phe Lys Ala Asp Gly Gly Pro Cys Lys Gly
 1                5                10                15
Lys Phe Pro Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
                20                25                30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
                35                40                45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50                55

<210> SEQ ID NO 45
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

-continued

<400> SEQUENCE: 45

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Lys Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Ala
 20 25 30

Phe Val Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 46

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 46

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Lys Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Pro
 20 25 30

Phe Val Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 47

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 47

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Lys Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Asn Thr
 20 25 30

Phe Val Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 48

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 48

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Gly
 1 5 10 15

Lys Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Gln Gly
 20 25 30

Phe Val Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

-continued

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 49
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 49

Met His Ser Phe Cys Ala Phe Lys Ala Glu Val Gly Pro Cys Arg Ala
1 5 10 15

Lys Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys His Leu
20 25 30

Phe Thr Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 50
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 50

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Gly
1 5 10 15

Lys Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Ala Gln
20 25 30

Phe Val Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 51
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 51

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Gly
1 5 10 15

Lys Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Ser
20 25 30

Phe His Tyr Gly Gly Cys Lys Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 52
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:

-continued

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 52

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ala Gly Pro Cys Arg Ala
 1 5 10 15
 Lys Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Ala
 20 25 30
 Phe Leu Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45
 Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 53

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 53

Met His Ser Phe Cys Ala Phe Lys Ala Asp Val Gly Pro Cys Arg Ala
 1 5 10 15
 Lys Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Ala
 20 25 30
 Phe Leu Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45
 Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 54

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 54

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ala Gly Pro Cys Arg Ala
 1 5 10 15
 Lys Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Thr Ala
 20 25 30
 Phe Phe Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45
 Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 55

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 55

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Pro Cys Arg Ala
 1 5 10 15
 Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Thr Arg
 20 25 30

-continued

Phe Pro Tyr Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 56
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 56

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Lys Ile Pro Arg Leu Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Pro
 20 25 30

Phe Ile Trp Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 57
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 57

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ala Gly Pro Cys Arg Ala
 1 5 10 15

Lys Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 58
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 58

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Lys Gly
 1 5 10 15

Ser Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Asn Val
 20 25 30

Phe Arg Tyr Gly Gly Cys Arg Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 59
 <211> LENGTH: 58
 <212> TYPE: PRT

-continued

```

<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 59
Met His Ser Phe Cys Ala Phe Lys Ala Asp Ala Gly Pro Cys Arg Ala
 1           5           10           15
Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Asp Thr
 20           25           30
Phe Leu Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50           55

```

```

<210> SEQ ID NO 60
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 60
Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Pro Cys Lys Gly
 1           5           10           15
Arg Phe Gly Arg Leu Phe Phe Asn Ile Phe Thr Arg Gln Cys Thr Ala
 20           25           30
Phe Asp Trp Gly Gly Cys Gly Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50           55

```

```

<210> SEQ ID NO 61
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 61
Met His Ser Phe Cys Ala Phe Lys Ala Asp Thr Gly Pro Cys Arg Ala
 1           5           10           15
Arg Phe Asp Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Ala
 20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50           55

```

```

<210> SEQ ID NO 62
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 62
Met His Ser Phe Cys Ala Phe Lys Ala Asp Thr Gly Pro Cys Arg Ala
 1           5           10           15
Arg Phe Asp Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu

```

-continued

20	25	30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu		
35	40	45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp		
50	55	

<210> SEQ ID NO 63
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 63

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ala Gly Pro Cys Arg Ala		
1	5	10
15		
Arg Phe Asp Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu		
20	25	30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu		
35	40	45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp		
50	55	

<210> SEQ ID NO 64
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Homo sapiens

<400> SEQUENCE: 64

Lys Pro Asp Phe Cys Phe Leu Glu Glu Asp Pro Gly Ile Cys Arg Gly		
1	5	10
15		
Tyr Ile Thr Arg Tyr Phe Tyr Asn Asn Gln Thr Lys Gln Cys Glu Arg		
20	25	30
Phe Lys Tyr Gly Gly Cys Leu Gly Asn Met Asn Asn Phe Glu Thr Leu		
35	40	45
Glu Glu Cys Lys Asn Ile Cys Glu Asp Gly		
50	55	

<210> SEQ ID NO 65
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 65

Lys Pro Asp Phe Cys Phe Leu Glu Glu Asp Thr Gly Pro Cys Arg Gly		
1	5	10
15		
Arg Phe Asp Arg Tyr Phe Tyr Asn Asn Gln Thr Lys Gln Cys Glu Thr		
20	25	30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Met Asn Asn Phe Glu Thr Leu		
35	40	45
Glu Glu Cys Lys Asn Ile Cys Glu Asp Gly		
50	55	

<210> SEQ ID NO 66
 <211> LENGTH: 58
 <212> TYPE: PRT

-continued

<213> ORGANISM: Homo sapiens

<400> SEQUENCE: 66

Gly Pro Ser Trp Cys Leu Thr Pro Ala Asp Arg Gly Leu Cys Arg Ala
 1 5 10 15
 Asn Glu Asn Arg Phe Tyr Tyr Asn Ser Val Ile Gly Lys Cys Arg Pro
 20 25 30
 Phe Lys Tyr Ser Gly Cys Gly Gly Asn Glu Asn Asn Phe Thr Ser Lys
 35 40 45
 Gln Glu Cys Leu Arg Ala Cys Lys Lys Gly
 50 55

<210> SEQ ID NO 67

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 67

Gly Pro Ser Trp Cys Leu Thr Pro Ala Asp Thr Gly Pro Cys Arg Ala
 1 5 10 15
 Arg Phe Asp Arg Phe Tyr Tyr Asn Ser Val Ile Gly Lys Cys Glu Pro
 20 25 30
 Phe Ile Tyr Gly Gly Cys Gly Gly Asn Glu Asn Asn Phe Thr Ser Lys
 35 40 45
 Gln Glu Cys Leu Arg Ala Cys Lys Lys Gly
 50 55

<210> SEQ ID NO 68

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Homo sapiens

<400> SEQUENCE: 68

Glu Thr Asp Ile Cys Lys Leu Pro Lys Asp Glu Gly Thr Cys Arg Asp
 1 5 10 15
 Phe Ile Leu Lys Trp Tyr Tyr Asp Pro Asn Thr Lys Ser Cys Ala Arg
 20 25 30
 Phe Trp Tyr Gly Gly Cys Gly Gly Asn Glu Asn Lys Phe Gly Ser Gln
 35 40 45
 Lys Glu Cys Glu Lys Val Cys Ala Pro Val
 50 55

<210> SEQ ID NO 69

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 69

Glu Thr Asp Ile Cys Lys Leu Pro Lys Asp Thr Gly Pro Cys Arg Ala
 1 5 10 15
 Arg Phe Asp Lys Trp Tyr Tyr Asp Pro Asn Thr Lys Ser Cys Glu Pro
 20 25 30
 Phe Val Tyr Gly Gly Cys Gly Gly Asn Glu Asn Lys Phe Gly Ser Gln
 35 40 45

-continued

Lys Glu Cys Glu Lys Val Cys Ala Pro Val
50 55

<210> SEQ ID NO 70
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Homo sapiens

<400> SEQUENCE: 70

Asn Ala Glu Ile Cys Leu Leu Pro Leu Asp Tyr Gly Pro Cys Arg Ala
1 5 10 15
Leu Leu Leu Arg Tyr Tyr Tyr Asp Arg Tyr Thr Gln Ser Cys Arg Gln
20 25 30
Phe Leu Tyr Gly Gly Cys Glu Gly Asn Ala Asn Asn Phe Tyr Thr Trp
35 40 45
Glu Ala Cys Asp Asp Ala Cys Trp Arg Ile
50 55

<210> SEQ ID NO 71
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 71

Asn Ala Glu Ile Cys Leu Leu Pro Leu Asp Thr Gly Pro Cys Arg Ala
1 5 10 15
Arg Phe Asp Arg Tyr Tyr Tyr Asp Arg Tyr Thr Gln Ser Cys Glu Gln
20 25 30
Phe Leu Tyr Gly Gly Cys Glu Gly Asn Ala Asn Asn Phe Tyr Thr Trp
35 40 45
Glu Ala Cys Asp Asp Ala Cys Trp Arg Ile
50 55

<210> SEQ ID NO 72
<211> LENGTH: 61
<212> TYPE: PRT
<213> ORGANISM: Homo sapiens

<400> SEQUENCE: 72

Val Pro Lys Val Cys Arg Leu Gln Val Ser Val Asp Asp Gln Cys Glu
1 5 10 15
Gly Ser Thr Glu Lys Tyr Phe Phe Asn Leu Ser Ser Met Thr Cys Glu
20 25 30
Lys Phe Phe Ser Gly Gly Cys His Arg Asn Arg Ile Glu Asn Arg Phe
35 40 45
Pro Asp Glu Ala Thr Cys Met Gly Phe Cys Ala Pro Lys
50 55 60

<210> SEQ ID NO 73
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 73

-continued

```

Val Pro Lys Val Cys Arg Leu Gln Val Glu Thr Gly Pro Cys Arg Gly
 1          5          10          15
Lys Phe Glu Lys Tyr Phe Phe Asn Leu Ser Ser Met Thr Cys Glu Thr
          20          25          30
Phe Val Tyr Gly Gly Cys Glu Gly Asn Arg Asn Arg Phe Pro Asp Glu
          35          40          45
Ala Thr Cys Met Gly Phe Cys Ala Pro Lys
          50          55

```

```

<210> SEQ ID NO 74
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Homo sapiens

```

```

<400> SEQUENCE: 74

```

```

Ile Pro Ser Phe Cys Tyr Ser Pro Lys Asp Glu Gly Leu Cys Ser Ala
 1          5          10          15
Asn Val Thr Arg Tyr Tyr Phe Asn Pro Arg Tyr Arg Thr Cys Asp Ala
          20          25          30
Phe Thr Tyr Thr Gly Cys Gly Gly Asn Asp Asn Asn Phe Val Ser Arg
          35          40          45
Glu Asp Cys Lys Arg Ala Cys Ala Lys Ala
          50          55

```

```

<210> SEQ ID NO 75
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 75

```

```

Ile Pro Ser Phe Cys Tyr Ser Pro Lys Asp Thr Gly Pro Cys Arg Ala
 1          5          10          15
Arg Phe Thr Arg Tyr Tyr Phe Asn Pro Arg Tyr Arg Thr Cys Asp Ala
          20          25          30
Phe Thr Tyr Gly Gly Cys Gly Gly Asn Asp Asn Asn Phe Val Ser Arg
          35          40          45
Glu Asp Cys Lys Arg Ala Cys Ala Lys Ala
          50          55

```

```

<210> SEQ ID NO 76
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Homo sapiens

```

```

<400> SEQUENCE: 76

```

```

Lys Glu Asp Ser Cys Gln Leu Gly Tyr Ser Ala Gly Pro Cys Met Gly
 1          5          10          15
Met Thr Ser Arg Tyr Phe Tyr Asn Gly Thr Ser Met Ala Cys Glu Thr
          20          25          30
Phe Gln Tyr Gly Gly Cys Met Gly Asn Gly Asn Asn Phe Val Thr Glu
          35          40          45
Lys Glu Cys Leu Gln Thr Cys Arg Thr Val
          50          55

```

```

<210> SEQ ID NO 77

```

-continued

<211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 77

Lys Glu Asp Ser Cys Gln Leu Gly Tyr Glu Ala Gly Pro Cys Arg Gly
 1 5 10 15

Lys Phe Ser Arg Tyr Phe Tyr Asn Gly Thr Ser Met Ala Cys Glu Thr
 20 25 30

Phe Val Tyr Gly Gly Cys Gly Gly Asn Gly Asn Asn Phe Val Thr Glu
 35 40 45

Lys Glu Cys Leu Gln Thr Cys Arg Thr Val
 50 55

<210> SEQ ID NO 78
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Homo sapiens

<400> SEQUENCE: 78

Thr Val Ala Ala Cys Asn Leu Pro Ile Val Arg Gly Pro Cys Arg Ala
 1 5 10 15

Phe Ile Gln Leu Trp Ala Phe Asp Ala Val Lys Gly Lys Cys Val Leu
 20 25 30

Phe Pro Tyr Gly Gly Cys Gln Gly Asn Gly Asn Lys Phe Tyr Ser Glu
 35 40 45

Lys Glu Cys Arg Glu Tyr Cys Gly Val Pro
 50 55

<210> SEQ ID NO 79
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 79

Thr Val Ala Ala Cys Asn Leu Pro Ile Asp Thr Gly Pro Cys Arg Ala
 1 5 10 15

Arg Phe Gln Leu Trp Ala Phe Asp Ala Val Lys Gly Lys Cys Val Leu
 20 25 30

Phe Val Tyr Gly Gly Cys Gln Gly Asn Gly Asn Lys Phe Tyr Ser Glu
 35 40 45

Lys Glu Cys Arg Glu Tyr Cys Gly Val Pro
 50 55

<210> SEQ ID NO 80
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Homo sapiens

<400> SEQUENCE: 80

Val Arg Glu Val Cys Ser Glu Gln Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Met Ile Ser Arg Trp Tyr Phe Asp Val Thr Glu Gly Lys Cys Ala Pro
 20 25 30

-continued

Phe Phe Tyr Gly Gly Cys Gly Gly Asn Arg Asn Asn Phe Asp Thr Glu
 35 40 45

Glu Tyr Cys Met Ala Val Cys Gly Ser Ala
 50 55

<210> SEQ ID NO 81
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 81

Val Arg Glu Val Cys Ser Glu Gln Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Arg Phe Ser Arg Trp Tyr Phe Asp Val Thr Glu Gly Lys Cys Ala Pro
 20 25 30

Phe Phe Tyr Gly Gly Cys Gly Gly Asn Arg Asn Asn Phe Asp Thr Glu
 35 40 45

Glu Tyr Cys Met Ala Val Cys Gly Ser Ala
 50 55

<210> SEQ ID NO 82
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 82

Val Arg Glu Val Cys Ser Glu Gln Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Arg Phe Ser Arg Trp Tyr Phe Asp Val Thr Glu Gly Lys Cys Glu Pro
 20 25 30

Phe Ile Tyr Gly Gly Cys Gly Gly Asn Arg Asn Asn Phe Asp Thr Glu
 35 40 45

Glu Tyr Cys Met Ala Val Cys Gly Ser Ala
 50 55

<210> SEQ ID NO 83
 <211> LENGTH: 97
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated oligonucleotide
 <220> FEATURE:
 <221> NAME/KEY: misc_feature
 <222> LOCATION: 38, 43, 56, 61, 62
 <223> OTHER INFORMATION: n = a, t, g or c
 <220> FEATURE:
 <221> NAME/KEY: CDS
 <222> LOCATION: (7)...(90)

<400> SEQUENCE: 83

cctcct atg cat tcc ttc tgc gcc ttc aag gct ras rnt ggt nct tgt 48
 Met His Ser Phe Cys Ala Phe Lys Ala Xaa Xaa Gly Xaa Cys
 1 5 10

ara gct ans wtc nns cgt tks ttc ttc aac atc ttc acg cgt 90
 Xaa Xaa Xaa Xaa Xaa Arg Xaa Phe Phe Asn Ile Phe Thr Arg
 15 20 25

tccctcc 97

-continued

<210> SEQ ID NO 84
 <211> LENGTH: 26
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated DNA fragment

<400> SEQUENCE: 84

ggaggggaacg cgtgaagatg ttgaag

26

<210> SEQ ID NO 85
 <211> LENGTH: 28
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 10
 <223> OTHER INFORMATION: Xaa = Asp or Glu; Asn or Lys
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 11
 <223> OTHER INFORMATION: Xaa = Asp or Val; Ala or Gly; Ile or Thr; Asn
 or Ser
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 13
 <223> OTHER INFORMATION: Ser or Thr; Ala or Pro
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 15
 <223> OTHER INFORMATION: Xaa = Lys or Arg
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 16
 <223> OTHER INFORMATION: Xaa = Ala or Gly
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 17
 <223> OTHER INFORMATION: Xaa = Arg or Ser; Ile or Met; Asn or Lys; Thr
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 18
 <223> OTHER INFORMATION: Xaa = Phe or Ile
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 19
 <223> OTHER INFORMATION: Xaa = Glu or Gly; Asp or Ala; Val or Arg; Ser
 or Lys; Asn or Thr; Met or Ile; Gln or His; Leu or
 Pro; Phe or Tyr; Cys or Trp
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 21
 <223> OTHER INFORMATION: Xaa = leu or Trp; Phe or Cys

<400> SEQUENCE: 85

Met His Ser Phe Cys Ala Phe Lys Ala Xaa Xaa Gly Xaa Cys Xaa Xaa
 1 5 10 15

Xaa Xaa Xaa Arg Xaa Phe Phe Asn Ile Phe Thr Arg
 20 25

<210> SEQ ID NO 86
 <211> LENGTH: 77
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated oligonucleotide
 <220> FEATURE:

-continued

```

<221> NAME/KEY: CDS
<222> LOCATION: (7)...(54)
<220> FEATURE:
<221> NAME/KEY: misc_feature
<222> LOCATION: 19, 20, 22, 23, 28, 29
<223> OTHER INFORMATION: n = a, c, t, or g

<400> SEQUENCE: 86

cctcct acg cgt cag tgc nns nns ttc nnt trs ggt ggt tgt rrg ggt      48
  Thr Arg Gln Cys Xaa Xaa Phe Xaa Xaa Gly Gly Cys Xaa Gly
    1             5             10

aac cag gtcgtgctct ttagcacgac ctg      77
Asn Gln
  15

<210> SEQ ID NO 87
<211> LENGTH: 16
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 5, 6
<223> OTHER INFORMATION: Xaa = Glu or Gln; Asp or Gly; Val or Ala; Asn
or Lys; Thr or Met; Trp or Ile; His or Arg; Leu or
Pro; Tyr or Cys; Phe or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 8
<223> OTHER INFORMATION: Xaa = Phe or Leu; Ser or Pro; Ile or Thr; Ala
or Val; Tyr or Cys; His or Arg; Asn or Asp; Gly
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 9
<223> OTHER INFORMATION: Xaa = Tyr or Trp; Cys
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 13
<223> OTHER INFORMATION: Xaa = Lys or Arg; Glu or Gly

<400> SEQUENCE: 87

Thr Arg Gln Cys Xaa Xaa Phe Xaa Xaa Gly Gly Cys Xaa Gly Asn Gln
  1             5             10             15

<210> SEQ ID NO 88
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 88

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
  1             5             10             15

Arg Phe Glu Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
    20             25             30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
    35             40             45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
    50             55

<210> SEQ ID NO 89
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence

```

-continued

```

<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 89

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
 1             5             10             15
Arg Phe Gly Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20             25             30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 90
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 90

Met His Ser Phe Cys Ala Phe Lys Ala Glu Ser Gly Pro Cys Arg Ala
 1             5             10             15
Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20             25             30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 91
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 91

Met His Ser Phe Cys Ala Phe Lys Ala Asp Thr Gly Pro Cys Arg Gly
 1             5             10             15
Arg Phe Glu Arg Leu Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20             25             30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 92
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 92

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Ser Cys Arg Gly
 1             5             10             15
Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20             25             30

```

-continued

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 93
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 93

Met His Ser Phe Cys Ala Phe Lys Ala Glu Val Gly Pro Cys Arg Ala
 1 5 10 15

Ser Phe Pro Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 94
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 94

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
 1 5 10 15

Thr Phe Pro Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 95
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 95

Met His Ser Phe Cys Ala Phe Lys Ala Glu Val Gly Pro Cys Arg Ala
 1 5 10 15

Ser Phe His Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 96
 <211> LENGTH: 58

-continued

<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 96

Met His Ser Phe Cys Ala Phe Lys Ala Asp Thr Gly Pro Cys Arg Ala
1 5 10 15
Ser Phe Gly Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 97
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 97

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Gly
1 5 10 15
Met Phe Pro Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 98
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 98

Met His Ser Phe Cys Ala Phe Lys Ala Glu Gly Gly Pro Cys Arg Ala
1 5 10 15
Arg Phe Asn Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 99
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 99

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Arg Ala
1 5 10 15

-continued

Arg Ile Ser Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 100
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 100

Met His Ser Phe Cys Ala Phe Lys Ala Glu Gly Gly Pro Cys Arg Ala
 1 5 10 15

Lys Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 101
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 101

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Ala Cys Arg Gly
 1 5 10 15

Arg Phe Glu Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 102
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 102

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Pro Cys Arg Gly
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

-continued

```

<210> SEQ ID NO 103
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 103

Met His Ser Phe Cys Ala Phe Lys Ala Asp Thr Gly Pro Cys Arg Ala
 1             5             10             15

Ser Phe Pro Arg Leu Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20             25             30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 104
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 104

Met His Ser Phe Cys Ala Phe Lys Ala Glu Val Gly Pro Cys Arg Ala
 1             5             10             15

Arg Ile Gln Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20             25             30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 105
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 105

Met His Ser Phe Cys Ala Phe Lys Ala Glu Ser Gly Pro Cys Arg Ala
 1             5             10             15

Lys Phe Ala Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20             25             30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35             40             45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50             55

```

```

<210> SEQ ID NO 106
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 106

Met His Ser Phe Cys Ala Phe Lys Ala Glu Gly Gly Pro Cys Arg Ala

```

-continued

```

1           5           10           15
Lys Phe Ala Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
      20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
      50           55

```

```

<210> SEQ ID NO 107
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 107

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Thr Gly Ser Cys Arg Ala
  1           5           10           15
Lys Ile Glu Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
      20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
      50           55

```

```

<210> SEQ ID NO 108
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 108

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Pro Cys Lys Ala
  1           5           10           15
Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
      20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
      50           55

```

```

<210> SEQ ID NO 109
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 109

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Gly Gly Pro Cys Lys Gly
  1           5           10           15
Arg Phe Glu Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
      20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
      35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
      50           55

```


-continued

<210> SEQ ID NO 110
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 110

Met His Ser Phe Cys Ala Phe Lys Ala Glu Val Gly Ala Cys Lys Gly
1 5 10 15
Arg Phe His Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 111
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 111

Met His Ser Phe Cys Ala Phe Lys Ala Asp Gly Gly Pro Cys Arg Ala
1 5 10 15
Ser Phe Pro Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 112
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 112

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Ala Cys Arg Ala
1 5 10 15
Met Phe His Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 113
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 113

-continued

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Ala Cys Arg Ala
 1           5           10           15
Lys Phe Arg Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50           55

```

```

<210> SEQ ID NO 114
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 114

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Thr Cys Lys Ala
 1           5           10           15
Arg Phe Pro Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50           55

```

```

<210> SEQ ID NO 115
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 115

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Lys Gly
 1           5           10           15
Lys Ile Ala Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35           40           45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50           55

```

```

<210> SEQ ID NO 116
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 116

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Asp Ser Gly Ala Cys Lys Gly
 1           5           10           15
Lys Phe Glu Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20           25           30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35           40           45

```

-continued

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 117
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 117

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Ala Ala
1 5 10 15
Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 118
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 118

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Gly Ala
1 5 10 15
Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 119
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 119

Met His Ser Phe Cys Ala Phe Lys Ala Glu Thr Gly Pro Cys Asn Ala
1 5 10 15
Arg Phe Asp Arg Trp Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Ala
20 25 30
Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
35 40 45
Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50 55

<210> SEQ ID NO 120
<211> LENGTH: 96
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated DNA fragment

-continued

```

<220> FEATURE:
<221> NAME/KEY: misc_feature
<222> LOCATION: 37, 40, 41, 43, 44, 46, 47, 53, 62
<223> OTHER INFORMATION: n = a, t, c or g
<220> FEATURE:
<221> NAME/KEY: CDS
<222> LOCATION: (7)...(96)

<400> SEQUENCE: 120

cctcct atg cat tcc ttc tgc gcc ttc aag gct gas nct nnt nnt nnt      48
  Met His Ser Phe Cys Ala Phe Lys Ala Xaa Xaa Xaa Xaa Xaa
    1             5             10

ara rnt ara ttc gns crt tks ttc ttc aac atc ttc acg cgt cag tgc    96
Xaa Xaa Xaa Phe Xaa Xaa Xaa Phe Phe Asn Ile Phe Thr Arg Gln Cys
  15             20             25             30

<210> SEQ ID NO 121
<211> LENGTH: 71
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated DNA fragment
<220> FEATURE:
<221> NAME/KEY: misc_feature
<222> LOCATION: 19, 22, 23, 25, 26
<223> OTHER INFORMATION: n = a, t, c or g

<400> SEQUENCE: 121

cctcctccct ggttaccsny annannaccg taaacgaaag cctcgcactg acgcggaag    60

atggtgaaga a                                                         71

<210> SEQ ID NO 122
<211> LENGTH: 30
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 10
<223> OTHER INFORMATION: Xaa = Asp or Glu
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 11
<223> OTHER INFORMATION: Xaa = Thr or Ala; Ser or Pro
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 12
<223> OTHER INFORMATION: Xaa =Gly; Ala or Asp; Asn or Val; Ile or Thr;
  His or Arg; Cys or Pro; Phe or Ser; Tyr or Leu
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 13, 14
<223> OTHER INFORMATION: Xaa = Gly; Ala or Asp; Asn or Val; Ile or Thr;
  His or Arg; Leu or Pro; Tyr or Cys; Phe or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 15
<223> OTHER INFORMATION: Xaa = Lys or Arg
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 16
<223> OTHER INFORMATION: Xaa = Ala or Gly; Val or Asp; Ile or Thr; Asn
  or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 17
<223> OTHER INFORMATION: Xaa = Arg or Lys
<220> FEATURE:
<221> NAME/KEY: VARIANT

```

-continued

```

<222> LOCATION: 19
<223> OTHER INFORMATION: Xaa = Gly; Ala or Glu; Val or Asp
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 20
<223> OTHER INFORMATION: Xaa = Arg or Gln
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 21
<223> OTHER INFORMATION: Xaa = Leu or Trp; Phe or Cys

<400> SEQUENCE: 122

Met His Ser Phe Cys Ala Phe Lys Ala Xaa Xaa Xaa Xaa Xaa Xaa
 1           5           10           15

Xaa Phe Xaa Xaa Xaa Phe Phe Asn Ile Phe Thr Arg Gln Cys
      20           25           30

<210> SEQ ID NO 123
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 1, 2, 3, 4, 56, 57, 58
<223> OTHER INFORMATION: Xaa = any amino acid or may be absent
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 6, 7, 8, 9
<223> OTHER INFORMATION: Xaa = any amino acid
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 10
<223> OTHER INFORMATION: Xaa = Asp, Glu or Tyr
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 11
<223> OTHER INFORMATION: Xaa = Thr, Ala, Ser, Val or Asp
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 13
<223> OTHER INFORMATION: Xaa = Pro, Leu or Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 15
<223> OTHER INFORMATION: Xaa = Lys or Arg
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 16
<223> OTHER INFORMATION: Xaa = Ala or Gly
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 17
<223> OTHER INFORMATION: Xaa = Arg, Lys or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 18
<223> OTHER INFORMATION: Xaa = Phe or Ile
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 19
<223> OTHER INFORMATION: Xaa = Glu, Asp, Pro, Gly, Ser or Ile
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 21
<223> OTHER INFORMATION: Xaa = Phe, Trp or Tyr
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 22, 23
<223> OTHER INFORMATION: Xaa = Tyr or Phe
<220> FEATURE:
<221> NAME/KEY: VARIANT

```

-continued

```

<222> LOCATION: 24, 25, 26, 27, 28, 29, 41, 42, 44, 46, 47, 48, 49, 50,
52, 53, 54
<223> OTHER INFORMATION: Xaa = any amino acid
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 31
<223> OTHER INFORMATION: Xaa = Asp, Glu, Thr, Val, Gln or Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 32
<223> OTHER INFORMATION: Xaa = Thr, Ala, Glu, Pro or Gln
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 34
<223> OTHER INFORMATION: Xaa = Val, Ile, Thr, Leu, Phe, Tyr, His, Asp,
Ala or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 35
<223> OTHER INFORMATION: Xaa = Tyr or Trp
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 36
<223> OTHER INFORMATION: Xaa = Gly or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 39
<223> OTHER INFORMATION: Xaa = Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys or
Met
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 40
<223> OTHER INFORMATION: Xaa = Gly or Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 43
<223> OTHER INFORMATION: Xaa = Asn or Gly
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 45
<223> OTHER INFORMATION: Xaa = Phe or Tyr

<400> SEQUENCE: 123

Xaa Xaa Xaa Xaa Cys Xaa Xaa Xaa Xaa Xaa Xaa Gly Xaa Cys Xaa Xaa
 1             5             10             15

Xaa Xaa Xaa Arg Xaa Xaa Xaa Xaa Xaa Xaa Xaa Xaa Cys Xaa Xaa
 20             25             30

Phe Xaa Xaa Xaa Gly Cys Xaa Xaa Xaa Xaa Xaa Xaa Xaa Xaa Xaa Xaa
 35             40             45

Xaa Xaa Cys Xaa Xaa Xaa Cys Xaa Xaa Xaa
 50             55

<210> SEQ ID NO 124
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 124

Leu Pro Asn Val Cys Ala Phe Pro Met Glu Lys Gly Pro Cys Gln Thr
 1             5             10             15

Tyr Met Thr Arg Trp Phe Phe Asn Phe Glu Thr Gly Glu Cys Glu Leu
 20             25             30

Phe Ala Tyr Gly Gly Cys Gly Gly Asn Ser Asn Asn Phe Leu Arg Lys
 35             40             45

Glu Lys Cys Glu Lys Phe Cys Lys Phe Thr

```


-continued

<400> SEQUENCE: 128

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Asp Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 129

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 129

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Lys Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Glu
 20 25 30

Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 130

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 130

Thr Val Ala Ala Cys Asn Leu Pro Ile Val Thr Gly Pro Cys Arg Ala
 1 5 10 15

Arg Phe Gln Leu Trp Ala Phe Asp Ala Val Lys Gly Lys Cys Val Leu
 20 25 30

Phe Pro Tyr Gly Gly Cys Gln Gly Asn Gly Asn Lys Phe Tyr Ser Glu
 35 40 45

Lys Glu Cys Arg Glu Tyr Cys Gly Val Pro
 50 55

<210> SEQ ID NO 131

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 131

Thr Val Ala Ala Cys Asn Leu Pro Ile Val Thr Gly Pro Cys Arg Ala
 1 5 10 15

Arg Phe Gln Arg Trp Ala Phe Asp Ala Val Lys Gly Lys Cys Val Leu
 20 25 30

Phe Pro Tyr Gly Gly Cys Gln Gly Asn Gly Asn Lys Phe Tyr Ser Glu
 35 40 45

-continued

Lys Glu Cys Arg Glu Tyr Cys Gly Val Pro
50 55

<210> SEQ ID NO 132
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 132

Thr Val Ala Ala Cys Asn Leu Pro Ile Val Thr Gly Pro Cys Arg Ala
1 5 10 15

Arg Phe Gln Arg Trp Ala Phe Asp Ala Val Lys Gly Lys Cys Val Leu
20 25 30

Phe Val Tyr Gly Gly Cys Gln Gly Asn Gly Asn Lys Phe Tyr Ser Glu
35 40 45

Lys Glu Cys Arg Glu Tyr Cys Gly Val Pro
50 55

<210> SEQ ID NO 133
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 133

Thr Val Ala Ala Cys Asn Leu Pro Ile Glu Thr Gly Pro Cys Arg Ala
1 5 10 15

Arg Phe Asp Arg Trp Ala Phe Asp Ala Val Lys Gly Lys Cys Glu Thr
20 25 30

Phe Val Tyr Gly Gly Cys Gly Gly Asn Gly Asn Lys Phe Tyr Ser Glu
35 40 45

Lys Glu Cys Arg Glu Tyr Cys Gly Val Pro
50 55

<210> SEQ ID NO 134
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 134

Arg Pro Asp Phe Cys Leu Glu Pro Pro Tyr Thr Gly Pro Cys Lys Ala
1 5 10 15

Arg Phe Ile Arg Tyr Phe Tyr Asn Ala Lys Ala Gly Leu Cys Gln Thr
20 25 30

Phe Val Tyr Gly Gly Cys Arg Ala Lys Arg Asn Asn Phe Lys Ser Ala
35 40 45

Glu Asp Cys Met Arg Thr Cys Gly Gly Ala
50 55

<210> SEQ ID NO 135
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:

-continued

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 135

```

Arg Pro Asp Phe Cys Leu Glu Pro Pro Tyr Thr Gly Pro Cys Arg Ala
 1           5           10           15
Arg Phe Ile Arg Tyr Phe Tyr Asn Ala Lys Ala Gly Leu Cys Gln Thr
 20           25           30
Phe Val Tyr Gly Gly Cys Arg Ala Lys Arg Asn Asn Phe Lys Ser Ala
 35           40           45
Glu Asp Cys Met Arg Thr Cys Gly Gly Ala
 50           55

```

<210> SEQ ID NO 136

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 136

```

Arg Pro Asp Phe Cys Leu Glu Pro Pro Tyr Thr Gly Pro Cys Arg Ala
 1           5           10           15
Arg Phe Ile Arg Tyr Phe Tyr Asn Ala Lys Ala Gly Leu Cys Gln Thr
 20           25           30
Phe Val Tyr Gly Gly Cys Gly Ala Lys Arg Asn Asn Phe Lys Ser Ala
 35           40           45
Glu Asp Cys Met Arg Thr Cys Gly Gly Ala
 50           55

```

<210> SEQ ID NO 137

<211> LENGTH: 58

<212> TYPE: PRT

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 137

```

Arg Pro Asp Phe Cys Leu Glu Pro Pro Asp Thr Gly Pro Cys Arg Ala
 1           5           10           15
Arg Phe Asp Arg Tyr Phe Tyr Asn Ala Lys Ala Gly Leu Cys Glu Thr
 20           25           30
Phe Val Tyr Gly Gly Cys Gly Ala Lys Arg Asn Asn Phe Lys Ser Ala
 35           40           45
Glu Asp Cys Met Arg Thr Cys Gly Gly Ala
 50           55

```

<210> SEQ ID NO 138

<211> LENGTH: 13

<212> TYPE: DNA

<213> ORGANISM: Artificial Sequence

<220> FEATURE:

<223> OTHER INFORMATION: Synthetically generated DNA fragment

<220> FEATURE:

<221> NAME/KEY: CDS

<222> LOCATION: (1)...(12)

<400> SEQUENCE: 138

```

aaa agg cct cga g
Lys Arg Pro Arg
 1

```

-continued

```

<210> SEQ ID NO 139
<211> LENGTH: 4
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 139

Lys Arg Pro Arg
1

<210> SEQ ID NO 140
<211> LENGTH: 42
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated DNA fragment
<220> FEATURE:
<221> NAME/KEY: CDS
<222> LOCATION: (1)...(42)

<400> SEQUENCE: 140

aaa agg gaa gcg gcc gag cca tgg ggc gcc taa tag ctc gag      42
Lys Arg Glu Ala Ala Glu Pro Trp Gly Ala * * Leu Glu
1           5           10

<210> SEQ ID NO 141
<211> LENGTH: 12
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 141

Lys Arg Glu Ala Ala Glu Pro Trp Gly Ala Leu Glu
1           5           10

<210> SEQ ID NO 142
<211> LENGTH: 210
<212> TYPE: DNA
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated oligonucleotide
<220> FEATURE:
<221> NAME/KEY: CDS
<222> LOCATION: (1)...(198)

<400> SEQUENCE: 142

aaa agg gaa gcg gcc gag atg cat tcc ttc tgc gct ttc aaa gct gat      48
Lys Arg Glu Ala Ala Glu Met His Ser Phe Cys Ala Phe Lys Ala Asp
1           5           10           15

gac ggt cgg tgt aaa gct atc atg aaa cgt ttc ttc ttc aac att ttc      96
Asp Gly Pro Cys Lys Ala Ile Met Lys Arg Phe Phe Phe Asn Ile Phe
20           25           30

acg cgt cag tgc gag gaa ttc att tac ggt ggt tgt gaa ggt aac cag      144
Thr Arg Gln Cys Glu Glu Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln
35           40           45

aac cgg ttc gaa tct cta gag gaa tgt aag aag atg tgc act cgt gac      192
Asn Arg Phe Glu Ser Leu Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
50           55           60

ggc gcc taatagctcg ag      210
Gly Ala
65

```

-continued

<210> SEQ ID NO 143
 <211> LENGTH: 66
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide

 <400> SEQUENCE: 143

 Lys Arg Glu Ala Ala Glu Met His Ser Phe Cys Ala Phe Lys Ala Asp
 1 5 10 15

 Asp Gly Pro Cys Lys Ala Ile Met Lys Arg Phe Phe Phe Asn Ile Phe
 20 25 30

 Thr Arg Gln Cys Glu Glu Phe Ile Tyr Gly Gly Cys Glu Gly Asn Gln
 35 40 45

 Asn Arg Phe Glu Ser Leu Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50 55 60

 Gly Ala
 65

<210> SEQ ID NO 144
 <211> LENGTH: 132
 <212> TYPE: DNA
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated DNA fragment
 <220> FEATURE:
 <221> NAME/KEY: CDS
 <222> LOCATION: (7)...(132)
 <220> FEATURE:
 <221> NAME/KEY: misc_feature
 <222> LOCATION: 37, 40, 41, 43, 44, 46, 47, 53, 62, 115, 116, 118, 119,
 122
 <223> OTHER INFORMATION: n = a, t, c or g

 <400> SEQUENCE: 144

 cctcct atg cat tcc ttc tgc gcc ttc aag gct gas nct nnt nnt nnt 48
 Met His Ser Phe Cys Ala Phe Lys Ala Xaa Xaa Xaa Xaa Xaa
 1 5 10

 ara nnt ara ttc gns crt tks ttc ttc aac atc ttc acg cgt cag tgc 96
 Xaa Xaa Xaa Phe Xaa Xaa Xaa Phe Phe Asn Ile Phe Thr Arg Gln Cys
 15 20 25 30

 gag gct ttc gtt tac ggt nnt nnt rns ggt aac cag 132
 Glu Ala Phe Val Tyr Gly Xaa Xaa Xaa Gly Asn Gln
 35 40

<210> SEQ ID NO 145
 <211> LENGTH: 42
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 10
 <223> OTHER INFORMATION: Xaa = Asp or Glu
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 11
 <223> OTHER INFORMATION: Xaa = Thr or Ala; Ser or Pro
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 12
 <223> OTHER INFORMATION: Xaa =Gly; Ala or Asp; Asn or Val; Ile or Thr;
 His or Arg; Cys or Pro; Phe or Ser; Tyr or Leu

-continued

```

<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 13, 14, 37
<223> OTHER INFORMATION: Xaa = Gly; Ala or Asp; Asn or Val; Ile or Thr;
His or Arg; Leu or Pro; Tyr or Cys; Phe or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 15
<223> OTHER INFORMATION: Xaa = Lys or Arg
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 16
<223> OTHER INFORMATION: Xaa = Ala or Gly; Val or Asp; Ile or Thr; Asn
or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 17
<223> OTHER INFORMATION: Xaa = Arg or Lys
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 19
<223> OTHER INFORMATION: Xaa = Gly; Ala or Glu; Val or Asp
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 20
<223> OTHER INFORMATION: Xaa = Arg or Gln
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 21
<223> OTHER INFORMATION: Xaa = Leu or Trp; Phe or Cys
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 38
<223> OTHER INFORMATION: Xaa = Cys; Asp or Gly; Val or Ala; Thr or Asn;
Arg or Ile; Pro or His; Tyr or Leu; Phe or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 39
<223> OTHER INFORMATION: Xaa = Gly or Asp; Ala or Glu; Arg or Val; Lys
or Ser; Thr or Asn; Ile or Met

```

```

<400> SEQUENCE: 145

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Xaa Xaa Xaa Xaa Xaa Xaa Xaa
 1           5                10                15

```

```

Xaa Phe Xaa Xaa Xaa Phe Phe Asn Ile Phe Thr Arg Gln Cys Glu Ala
      20                25                30

```

```

Phe Val Tyr Gly Xaa Xaa Xaa Gly Asn Gln
 35                40

```

```

<210> SEQ ID NO 146
<211> LENGTH: 12
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 146

```

```

Glu Thr Gly Pro Cys Arg Ala Lys Phe Asp Arg Trp
 1           5                10

```

```

<210> SEQ ID NO 147
<211> LENGTH: 10
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

```

```

<400> SEQUENCE: 147

```

```

Glu Ala Phe Val Tyr Gly Gly Cys Gly Gly

```

-continued

1 5 10

<210> SEQ ID NO 148
<211> LENGTH: 12
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 148

Asp Thr Gly Pro Cys Arg Ala Arg Phe Asp Arg Phe
1 5 10

<210> SEQ ID NO 149
<211> LENGTH: 10
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide

<400> SEQUENCE: 149

Glu Ala Phe Ile Tyr Gly Gly Cys Glu Gly
1 5 10

<210> SEQ ID NO 150
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 10
<223> OTHER INFORMATION: Xaa = Asp, Glu, and Tyr
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 11
<223> OTHER INFORMATION: Xaa = Thr, Ala, Ser, Val, and Asp;
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 13
<223> OTHER INFORMATION: Xaa = Pro, Leu, and Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 15
<223> OTHER INFORMATION: Xaa = Arg and Lys
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 16
<223> OTHER INFORMATION: Xaa = Ala and Gly
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 17
<223> OTHER INFORMATION: Xaa = Arg, Lys, and Ser;
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 18
<223> OTHER INFORMATION: Xaa = Phe and Ile
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 19
<223> OTHER INFORMATION: Xaa = Glu, Asp, Pro, Gly, Ser, or Ile
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 22, 23
<223> OTHER INFORMATION: Xaa = Tyr or Phe
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 31
<223> OTHER INFORMATION: Xaa = Asp, Glu, Thr, Val, Gln, or Ala
<220> FEATURE:

-continued

```

<221> NAME/KEY: VARIANT
<222> LOCATION: 32
<223> OTHER INFORMATION: Xaa = Thr, Ala, Glu, Pro, and Gln
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 34
<223> OTHER INFORMATION: Xaa = Val, Ile, Thr, Leu, Phe, Tyr, His, Asp,
Ala, or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 35
<223> OTHER INFORMATION: Xaa = Tyr or Trp
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 39
<223> OTHER INFORMATION: Xaa = Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys,
Or Met
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 40
<223> OTHER INFORMATION: Xaa = Gly and Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 43
<223> OTHER INFORMATION: Xaa = Asn and Gly
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 45
<223> OTHER INFORMATION: Xaa = Phe and Tyr

```

```

<400> SEQUENCE: 150

```

```

Met His Ser Phe Cys Ala Phe Lys Ala Xaa Xaa Gly Xaa Cys Xaa Xaa
  1             5             10             15

```

```

Xaa Xaa Xaa Arg Trp Xaa Xaa Asn Ile Phe Thr Arg Gln Cys Xaa Xaa
  20             25             30

```

```

Phe Xaa Xaa Gly Gly Cys Xaa Xaa Asn Gln Xaa Arg Xaa Glu Ser Leu
  35             40             45

```

```

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
  50             55

```

```

<210> SEQ ID NO 151
<211> LENGTH: 58
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 10
<223> OTHER INFORMATION: Xaa = Asp or Glu or Tyr
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 11
<223> OTHER INFORMATION: Xaa = Thr, Ala, Ser, Val, or Asp
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 13
<223> OTHER INFORMATION: Xaa = Pro, Leu or Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 15
<223> OTHER INFORMATION: Xaa = Arg or Lys
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 16
<223> OTHER INFORMATION: Ala or Gly
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 17
<223> OTHER INFORMATION: Xaa = Arg, Lys or Ser
<220> FEATURE:

```

-continued

```

<221> NAME/KEY: VARIANT
<222> LOCATION: 18
<223> OTHER INFORMATION: Xaa = Phe or Ile
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 19
<223> OTHER INFORMATION: Xaa = Glu, Asp, Pro, Gly Ser or Ile
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 21
<223> OTHER INFORMATION: Xaa = Phe, Trp or Tyr
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 22, 23
<223> OTHER INFORMATION: Xaa = Tyr or Phe
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 31, 32
<223> OTHER INFORMATION: Xaa = Glu, Gln; Gln, Gln; or Gln, Glu
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 34
<223> OTHER INFORMATION: Xaa = Val, Ile, Thr, Leu, Phe, Tyr, His, Asp,
Ala or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 35
<223> OTHER INFORMATION: Xaa = Tyr or Trp
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 39
<223> OTHER INFORMATION: Xaa = Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys,
Phe, Asn, His or Met
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 40
<223> OTHER INFORMATION: Xaa = Gly or Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 43
<223> OTHER INFORMATION: Xaa = Asn or Gly
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 45
<223> OTHER INFORMATION: Xaa = Phe or Tyr

<400> SEQUENCE: 151

Met His Ser Phe Cys Ala Phe Lys Ala Xaa Xaa Gly Xaa Cys Xaa Xaa
 1           5           10          15

Xaa Xaa Xaa Arg Xaa Xaa Xaa Asn Ile Phe Thr Arg Gln Cys Xaa Xaa
 20          25          30

Phe Xaa Xaa Gly Gly Cys Xaa Xaa Asn Gln Xaa Arg Xaa Glu Ser Leu
 35          40          45

Glu Glu Cys Lys Lys Met Cys Thr Arg Asp
 50          55

<210> SEQ ID NO 152
<211> LENGTH: 56
<212> TYPE: PRT
<213> ORGANISM: Artificial Sequence
<220> FEATURE:
<223> OTHER INFORMATION: Synthetically generated peptide
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 31, 32
<223> OTHER INFORMATION: Xaa = Glu, Gln; Gln, Gln; or Gln, Glu
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 34
<223> OTHER INFORMATION: Xaa = Val Ile, Thr, Leu, Phe, Tyr, His, Asp,
Ala, or Ser

```


-continued

<220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 39
 <223> OTHER INFORMATION: Xaa = Gly, Gly, Asp, Arg, Ala, Gln, Leu, Lys,
 Phe, Asn, His or Met

<400> SEQUENCE: 152

Met His Ser Phe Cys Ala Phe Lys Ala Asp Asp Gly Pro Cys Lys Ala
 1 5 10 15

Arg Phe Glu Arg Phe Phe Phe Asn Ile Phe Thr Arg Gln Cys Xaa Xaa
 20 25 30

Phe Xaa Tyr Gly Gly Cys Xaa Gly Asn Gln Asn Arg Phe Glu Ser Leu
 35 40 45

Glu Glu Lys Met Cys Thr Arg Asp
 50 55

<210> SEQ ID NO 153
 <211> LENGTH: 58
 <212> TYPE: PRT
 <213> ORGANISM: Artificial Sequence
 <220> FEATURE:
 <223> OTHER INFORMATION: Synthetically generated peptide
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 1, 2, 3, 4, 56, 57, 58
 <223> OTHER INFORMATION: Xaa = any amino acid or may be absent
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 6, 7, 8, 9
 <223> OTHER INFORMATION: Xaa = any amino acid
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 10
 <223> OTHER INFORMATION: Xaa = Asp or Glu
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 11
 <223> OTHER INFORMATION: Xaa = any amino acid other than Asp or Glu
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 13
 <223> OTHER INFORMATION: Xaa = Pro, Leu or Ala
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 15
 <223> OTHER INFORMATION: Xaa = Arg or Lys
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 16
 <223> OTHER INFORMATION: Xaa = Ala or Gly
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 17
 <223> OTHER INFORMATION: Xaa = Arg, Lys or Ser
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 18
 <223> OTHER INFORMATION: Xaa = Phe or Ile
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 19
 <223> OTHER INFORMATION: Xaa = Glu, Gln, Asp, Pro, Gly, Ser or Ile
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 21
 <223> OTHER INFORMATION: Xaa = Phe, Trp or Tyr
 <220> FEATURE:
 <221> NAME/KEY: VARIANT
 <222> LOCATION: 22, 23
 <223> OTHER INFORMATION: Xaa = Tyr or Phe
 <220> FEATURE:

-continued

```

<221> NAME/KEY: VARIANT
<222> LOCATION: 24, 25, 26, 27, 28, 29, 36, 41, 42, 44, 46, 47, 48, 49,
50, 52, 53, 54,
<223> OTHER INFORMATION: Xaa = any amino acid
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 31
<223> OTHER INFORMATION: Xaa = Asp, Glu, Thr, Val, Gln or Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 32
<223> OTHER INFORMATION: Xaa = Thr, Ala, Glu, Pro or Gln
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 34
<223> OTHER INFORMATION: Xaa = Val, Ile, Thr, Leu, Phe, Tyr, His, Asp,
Ala or Ser
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 35
<223> OTHER INFORMATION: Xaa = Tyr or Trp
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 39
<223> OTHER INFORMATION: Xaa = Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys
or Met
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 40
<223> OTHER INFORMATION: Xaa = Gly or Ala
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 43
<223> OTHER INFORMATION: Xaa = Asn or Gly
<220> FEATURE:
<221> NAME/KEY: VARIANT
<222> LOCATION: 45
<223> OTHER INFORMATION: Xaa = Phe or Tyr

<400> SEQUENCE: 153

Xaa Xaa Xaa Xaa Cys Xaa Xaa Xaa Xaa Xaa Xaa Gly Xaa Cys Xaa Xaa
 1             5             10             15

Xaa Xaa Xaa Arg Xaa Xaa Xaa Xaa Xaa Xaa Xaa Xaa Cys Xaa Xaa
 20             25             30

Phe Xaa Xaa Xaa Gly Cys Xaa Xaa Xaa Xaa Xaa Xaa Xaa Xaa Xaa
 35             40             45

Xaa Xaa Cys Xaa Xaa Xaa Cys Xaa Xaa Xaa
 50             55

```

1.-11. (canceled)

12. A method of assaying a sample for plasmin comprising: contacting the sample with a polypeptide comprising a non-naturally occurring Kunitz domain that has the formula:

Met-His-Ser-Phe-Cys-Ala-Phe-Lys-Ala-Xaa10-
Xaa11-Gly-Xaa13-Cys-Xaa15-Xaa16-Xaa17-Xaa18-
Xaa19-Arg-Trp-Xaa22-Xaa23-Asn-Ile-Phe-Thr-Arg-
Gln-Cys-Xaa31-Xaa32-Phe-Xaa34-Xaa35-Gly-Gly-
Cys-Xaa39-Xaa40-Asn-Gln-Xaa43-Arg-Xaa45-Glu-
Ser-Leu-Glu-Glu-Cys-LysLys Met-Cys-Thr-Arg-Asp,

wherein Xaa10 is selected from the group consisting of Asp, Glu, and Tyr;

Xaa11 is selected from the group consisting of Thr, Ala, Ser, Val, and Asp;

Xaa13 is selected from the group consisting of Pro, Leu, and Ala;

Xaa15 is selected from the group consisting of Arg and Lys;

Xaa16 is selected from the group consisting of Ala and Gly;
Xaa17 is selected from the group consisting of Arg, Lys, and Ser;

Xaa18 is selected from the group consisting of Phe and Ile;
Xaa19 is selected from the group consisting of Glu, Asp, Pro, Gly, Ser, and Ile;

Xaa22 is selected from the group consisting of Tyr and Phe;
Xaa23 is selected from the group consisting of Tyr and Phe;
Xaa31 is selected from the group consisting of Asp, Glu, Thr, Val, Gln, and Ala;

Xaa32 is selected from the group consisting of Thr, Ala, Glu, Pro, and Gln;

Xaa34 is selected from the group consisting of Val, Ile, Thr, Leu, Phe, Tyr, His, Asp, Ala, and Ser;

Xaa35 is selected from the group consisting of Tyr and Trp;

Xaa39 is selected from the group consisting of Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys, and Met;

Xaa40 is selected from the group consisting of Gly and Ala; Xaa43 is selected from the group consisting of Asn and Gly; and

Xaa45 is selected from the group consisting of Phe and Tyr; and

detecting the presence of a complex of said polypeptide and plasmin in the sample.

13. The method according to claim **12**, wherein said polypeptide comprises a Kunitz domain having the formula:

Met-His-Ser-Phe-Cys-Ala-Phe-Lys-Ala-Xaa10-
Xaa11-Gly-Xaa13-Cys-Xaa15-Xaa16-Xaa17-Xaa18-
Xaa19-Arg-Trp-Xaa22-Xaa23-Asn-Ile-Phe-Thr-Arg-
Gln-Cys-Xaa31-Xaa32-Phe-Xaa34-Xaa35-Gly-Gly-
Cys-Xaa39-Xaa40-Asn-Gln-Xaa43-Arg-Xaa45-Glu-
Ser-Leu-Glu-Glu-Cys-Lys-Lys-Met-Cys-Thr-Arg-
Asp,

wherein Xaa10 is selected from the group consisting of Asp and Glu;

Xaa11 is selected from the group consisting of Thr, Ala, Ser, Val, and Asp;

Xaa13 is selected from the group consisting of Pro, Leu, and Ala;

Xaa15 is selected from the group consisting of Arg and Lys;

Xaa16 is selected from the group consisting of Ala and Gly; Xaa17 is selected from the group consisting of Arg, Lys, and Ser;

Xaa18 is selected from the group consisting of Phe and Ile; Xaa19 is selected from the group consisting of Glu, Asp, Pro, Gly, Ser, and Ile;

Xaa22 is Phe; Xaa23 is Phe; Xaa31 is selected from the group consisting of Asp, Glu, Thr, Val, Gln, and Ala;

Xaa32 is selected from the group consisting of Thr, Ala, Glu, Pro, and Gln;

Xaa34 is selected from the group consisting of Val, Ile, Thr, Leu, Phe, Tyr, His, Asp, Ala, and Ser;

Xaa35 is selected from the group consisting of Tyr and Trp; Xaa39 is selected from the group consisting of Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys, and Met;

Xaa40 is selected from the group consisting of Gly and Ala; Xaa43 is selected from the group consisting of Asn and Gly; and Xaa45 is selected from the group consisting of Phe and Tyr.

14. The method according to claim **12**, wherein said polypeptide comprises a Kunitz domain having a sequence selected from the group consisting of:

(SEQ ID NO: 40)
MHSFCAPKAETGPCRARFDRWFFNIFTRQCEEFIYGGCEGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 41)
MHSFCAPKAESGPCRARFDRWFFNIFTRQCEEFIYGGCEGNQNRFS
ESLBECKKMCTRD;

(SEQ ID NO: 43)
MHSFCAPKAEGGPCRAKFRWFFNIFTRQCEEFIYGGCEGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 45)
MHSFCAPKAETGPCRAKFRWFFNIFTRQCEAEVYGGCGGNQNRFSLEE
CKKMCTRD;

-continued

(SEQ ID NO: 46)
MHSFCAPKAETGPCRAKFRWFFNIFTRQCEPFVYGGCEGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 47)
MHSFCAPKAETGPCRAKFRWFFNIFTRQCNTPVYGGCGGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 48)
MHSFCAPKAETGPCRGKFRWFFNIFTRQCQPFVYGGCGGNQNRFSLEEC
KKMCTRD;

(SEQ ID NO: 49)
MHSFCAPKAIEVGPCRAKFRWFFNIFTRQCHLFTYGGCGGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 50)
MHSFCAPKAETGPCRGKFRWFFNIFTRQCAQPFVYGGCEGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 51)
MHSFCAPKAETGPCRGKFRWFFNIFTRQCESFHYGGCKGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 55)
MHSFCAPKADSGPCRARFDRWFFNIFTRQTRYPYGGCGGNQNRFSLEEC
KKMCTRD;

(SEQ ID NO: 58)
MHSFCAPKAETGPCKGSFDRWFFNIFTRQCNVFRYGGCGGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 117)
MHSFCAPKAETGPCAARFDRWFFNIFTRQCEEFIYGGCEGNQNRFSLEE
CKKMCTRD;

(SEQ ID NO: 118)
MHSFCAPKAETGPCGARFDRWFFNIFTRQCEEFIYGGCEGNQNPYSLEE
CKKMCTRD;

and

(SEQ ID NO: 119)
MHSFCAPKAETGPCNARFDRWFFNIFTRQCEAFIYGGCEGNQNRFSLEE
CKKMCTRD.

15. The method according to claim **12**, wherein said polypeptide comprises a Kunitz domain having the sequence:

(SEQ ID NO: 40)
MHSFCAPKAETGPCRARFDRWFFNIFTRQCEEFIYGGCEGNQNRFSLEE
CKKMCTRD.

16. A method of assaying a sample for plasmin comprising: contacting the sample with a polypeptide comprising a non-naturally occurring Kunitz domain that has the formula:

Met-His-Ser-Phe-Cys-Ala-Phe-Lys-Ala-Xaa10-
Xaa11-Gly-Xaa13-Cys-Xaa15-Xaa16-Xaa17-Xaa18-
Xaa19-Arg-Xaa21-Xaa22-Xaa23-Asn-Ile-Phe-Thr-
Arg-Gln-Cys-Xaa31-Xaa32-Phe-Xaa34-Xaa35-Gly-
Gly-Cys-Xaa39-Xaa40-Asn-Gln-Xaa43-Arg-Xaa45-
Glu-Ser-Leu-Glu-Glu-Cys-Lys-Lys-Met-Cys-Thr-
Arg-Asp,

wherein Xaa10 is selected from the group consisting of Asp, Glu, and Tyr;

Xaa11 is selected from the group consisting of Thr, Ala, Ser, Val, and Asp;

Xaa13 is selected from the group consisting of Pro, Leu, and Ala;

Xaa15 is selected from the group consisting of Arg and Lys;

Xaa16 is selected from the group consisting of Ala and Gly;
 Xaa17 is selected from the group consisting of Arg, Lys, and Ser;
 Xaa18 is selected from the group consisting of Phe and Ile;
 Xaa19 is selected from the group consisting of Glu, Asp, Pro, Gly, Ser, and Ile; Xaa21 is selected from the group consisting of Phe, Trp, and Tyr;
 Xaa22 is selected from the group consisting of Tyr and Phe;
 Xaa23 is selected from the group consisting of Tyr and Phe;
 Xaa31-Xaa32 is selected from the group consisting of Glu-Gln, Gln-Gln, and Gln-Glu;
 Xaa34 is selected from the group consisting of Val, Ile, Thr, Leu, Phe, Tyr, His, Asp, Ala, and Ser;
 Xaa35 is selected from the group consisting of Tyr and Trp;
 Xaa39 is selected from the group consisting of Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys, Phe, Asn, His and Met;
 Xaa40 is selected from the group consisting of Gly and Ala;
 Xaa43 is selected from the group consisting of Asn and Gly; and
 Xaa45 is selected from the group consisting of Phe and Tyr.
17. The method according to claim 16, wherein said polypeptide comprises a Kunitz domain having the formula:

Met-His-Ser-Phe-Cys-Ala-Phe-Lys-Ala-Asp-Asp-Gly-Pro-Cys-Lys-Ala-Arg-Phe-Glu-Arg-Phe-Phe-Phe-Asn-Ile-Phe-Thr-Arg-Gln-Cys-Xaa31-Xaa32-Phe-Xaa34-Tyr-Gly-Gly-Cys-Xaa39-Gly-Asn-Gln-Asn-Arg-Phe-Glu-Ser-Leu-Glu-Glu-Lys-Met-Cys-Thr-Arg-Asp,

wherein Xaa31-Xaa32 is selected from the group consisting of Glu-Gln, Gln-Gln, and Gln-Glu;
 Xaa34 is selected from the group consisting of Val, Ile, Thr, Leu, Phe, Tyr, His, Asp, Ala, and Ser; and
 Xaa39 is selected from the group consisting of Glu, Gly, Asp, Arg, Ala, Gln, Leu, Lys, Phe, Asn, His and Met.
18. The method according to claim 16, wherein said polypeptide comprises a Kunitz domain having a sequence selected from the group consisting of:

(SEQ ID NO: 3)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFTYGGCRGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 4)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFYGGCDGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 5)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQPHYGGCDGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 6)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFDYGGCAGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 7)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQEFYGGCDGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 8)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCQFYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 10)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCQFYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 11)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFTYGGCQGNQNRFSLEE
 CKKMCTRD;

-continued

(SEQ ID NO: 13)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFIYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 14)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFIYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 15)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFIYGGCFGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 16)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCQPHYGGCEGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 17)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFVYGGCAGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 18)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFLYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 19)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFIYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 20)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFVYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 23)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFTYGGCMGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 24)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFSYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 26)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFSYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 27)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFAYGGCAGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 28)
 MHSFCAFKADDGPCKAFERFFFNIFTRQCEQFIYGGCVGNQNRFSLEEC
 KLMCTRD;

(SEQ ID NO: 31)
 MHSFCAFKADDGPCKASFERFFFNIFTRQCEQFTYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 33)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFFYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 34)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFTYGGCQGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 35)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFTYGGCMGNQNRFSLEE
 CKKMCTRD;

(SEQ ID NO: 36)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFTYGGCQGNQNRFSLEE
 CKKMCTRD;

and

(SEQ ID NO: 37)
 MHSFCAFKADDGPCKARFERFFFNIFTRQCEQFVYGGCQGNQNRFSLEE
 CKKMCTRD.

19. The method according to claim 18, wherein said polypeptide comprises a Kunitz domain having the sequence:

(SEQ ID NO: 4)
 MHSFCAFKADDGPCKARFEREFFNIFTRQCEQFYGGCDGNQNRFSLEE
 CKKMCTRD.

20. The method according to claim 12, wherein said polypeptide has a K_i for human plasmin of 100 pM or less.

21. The method according to claim 16, wherein said polypeptide has a K_i for human plasmin of 100 pM or less.

22. The method according to claim 12, wherein said polypeptide further comprises one or more amino acids upstream of the Kunitz domain.

23. The method according to claim 22, wherein the one or more amino acids are involved with processing by a recombinant host cell.

24. The method according to claim 15, wherein said polypeptide further comprises one or more amino acids upstream of the Kunitz domain.

25. The method according to claim 24, wherein the one or more amino acids are involved with processing by a recombinant host cell.

26. The method according to claim 16, wherein said polypeptide further comprises one or more amino acids upstream of the Kunitz domain.

27. The method according to claim 26, wherein the one or more amino acids are involved with processing by a recombinant host cell.

28. The method according to claim 19, wherein said polypeptide further comprises one or more amino acids upstream of the Kunitz domain.

29. The method according to claim 12, wherein the polypeptide comprises a label

30. The method according to claim 29, wherein the label is selected from the group consisting of a radioisotope, a fluorophore, an enzyme, a co-enzyme, an enzyme substrate, an electron-dense compound and an agglutinable particle.

31. The method according to claim 16, wherein the polypeptide comprises a label

32. The method according to claim 31, wherein the label is selected from the group consisting of a radioisotope, a fluorophore, an enzyme, a co-enzyme, an enzyme substrate, an electron-dense compound and an agglutinable particle.

* * * * *