

(19) **DANMARK**



Patent- og
Varemærkestyrelsen

(10) **DK/EP 3666902 T3**

(12) Oversættelse af
europæisk patentskrift

-
- (51) Int.Cl.: **C 12 Q 1/6827 (2018.01)** **C 12 Q 1/6883 (2018.01)**
- (45) Oversættelsen bekendtgjort den: **2024-09-09**
- (80) Dato for Den Europæiske Patentmyndigheds bekendtgørelse om meddelelse af patentet: **2024-07-03**
- (86) Europæisk ansøgning nr.: **19216148.7**
- (86) Europæisk indleveringsdag: **2016-05-20**
- (87) Den europæiske ansøgnings publiceringsdag: **2020-06-17**
- (30) Prioritet: **2015-05-22 US 201562165593 P** **2015-12-04 US 201562263320 P**
- (62) Stamansøgningsnr: **16738215.9**
- (84) Designerede stater: **AL AT BE BG CH CY CZ DE DK EE ES FI FR GB GR HR HU IE IS IT LI LT LU LV MC MK MT NL NO PL PT RO RS SE SI SK SM TR**
- (73) Patenthaver: **MEDICOVER PUBLIC CO LTD, Neas Engomis 31 , Engomi, Nicosia 2409, Cypern**
- (72) Opfinder: **KOUMBARIS, George, 101, Makariou III Street, 2565 Lithrodontas, Cypern**
PATSALIS, Philippos, 34, Gregoriou Xenopoulou Engomi, 2402 Nicosia, Cypern
KYPRI, Elena, 4, Saphous Street, Flat 204 Strovolos, 2007 Nico, Cypern
TSANGARAS, Kyriakos, 7, Kleonos, 3085 Limassol, Cypern
ACHILLEOS, Acilleas, 6 Ikarou St. , Yermasoyia, 4045 Limassol, Cypern
PAPAGEORGIOU, Elisavet, A., Kykladon 80 , Flat 201, 2301 Nicosia, Cypern
MINA, Petros, 133, Athalassas Avenue , Dasoupolis, 2015 Nicosia, Cypern
- (74) Fuldmægtig i Danmark: **AWA Denmark A/S, Strandgade 56, 1401 København K, Danmark**
- (54) Benævnelse: **MULTIPLEX PARALLEL ANALYSE AF MÅLRETTEDE GENOMISKE REGIONER TIL IKKE-INVASIV PRÆNATAL TESTNING**
- (56) Fremdragne publikationer:
PARAYRE ET AL: "Easy DNA extraction method and optimisation of PCR-Temporal Temperature Gel Electrophoresis to identify the predominant high and low GC-content bacteria from dairy products", JOURNAL OF MICROBIOLOGICAL METHODS, ELSEVIER, AMSTERDAM, NL, vol. 69, no. 3, 5 May 2007 (2007-05-05), pages 431-441, XP022059909, ISSN: 0167-7012, DOI: 10.1016/J.MIMET.2007.02.011
Y. LI ET AL: "GC Content-Based Pan-Pox Universal PCR Assays for Poxvirus Detection", JOURNAL OF CLINICAL MICROBIOLOGY, vol. 48, no. 1, 11 November 2009 (2009-11-11), pages 268-276, XP055064098, ISSN: 0095-1137, DOI: 10.1128/JCM.01697-09
GARY J. W. LIAO ET AL: "Noninvasive Prenatal Diagnosis of Fetal Trisomy 21 by Allelic Ratio Analysis Using Targeted Massively Parallel Sequencing of Maternal Plasma DNA", PLOS ONE, vol. 7, no. 5, 29 May 2012 (2012-05-29), page e38154, XP055062111, DOI: 10.1371/journal.pone.0038154
ELISAVET A PAPAGEORGIOU ET AL: "Non-invasive prenatal diagnosis of aneuploidies: new technologies and clinical applications", GENOME MED, vol. 4, no. 5, 1 January 2012 (2012-01-01) , page 46, XP055305603, London, UK ISSN: 1756-994X, DOI: 10.1186/gm345
PARAYRE ET AL: "Easy DNA extraction method and optimisation of PCR-Temporal Temperature Gel

Fortsættes ...

Electrophoresis to identify the predominant high and low GC-content bacteria from dairy products", JOURNAL OF MICROBIOLOGICAL METHODS, ELSEVIER, AMSTERDAM, NL, vol. 69, no. 3, 5 May 2007 (2007-05-05), pages 431 - 441, XP022059909, ISSN: 0167-7012, DOI: 10.1016/J.MIMET.2007.02.011

Y. LI ET AL: "GC Content-Based Pan-Pox Universal PCR Assays for Poxvirus Detection", JOURNAL OF CLINICAL MICROBIOLOGY, vol. 48, no. 1, 11 November 2009 (2009-11-11), pages 268 - 276, XP055064098, ISSN: 0095-1137, DOI: 10.1128/JCM.01697-09

GARY J. W. LIAO ET AL: "Noninvasive Prenatal Diagnosis of Fetal Trisomy 21 by Allelic Ratio Analysis Using Targeted Massively Parallel Sequencing of Maternal Plasma DNA", PLOS ONE, vol. 7, no. 5, 29 May 2012 (2012-05-29), pages e38154, XP055062111, DOI: 10.1371/journal.pone.0038154

ELISAVET A PAPAGEORGIOU ET AL: "Non-invasive prenatal diagnosis of aneuploidies: new technologies and clinical applications", GENOME MED, vol. 4, no. 5, 1 January 2012 (2012-01-01), London, UK, pages 46, XP055305603, ISSN: 1756-994X, DOI: 10.1186/gm345

DESCRIPTION

Description

Technical field of the invention

[0001] The present invention is in the field of non-invasive prenatal testing, in particular in the field of testing the risk of a chromosomal abnormality in fetal DNA. The present invention relates to a kit for performing a method of testing the risk of a chromosomal abnormality in fetal DNA.

Background of the Invention

[0002] The discovery of free fetal DNA (ffDNA) in maternal circulation (Lo, Y.M. et al. (1997) Lancet 350:485-487) was a landmark towards the development of non-invasive prenatal testing for fetal aneuploidies and has opened up new possibilities in the clinical setting. ffDNA has been successfully used for the determination of fetal sex and fetal Rhesus D status in maternal plasma (see e.g., Bianchi, D. et al. (2005) Obstet. Gynecol. 106:841-844; Lo, Y.M. et al. (1998) N. Engl. J. Med. 339:1734-1738; US Patent No. 6,258,540; PCT Publication WO 91/07660). These methods have become routine tests in a number of diagnostic laboratories worldwide. However, direct analysis of the limited amount of ffDNA in the presence of an excess of maternal DNA is a great challenge for Non-invasive Prenatal Testing (NIPT) assessment of fetal aneuploidies.

[0003] The percentage of ffDNA in the maternal circulation was originally estimated to be about 3-6% of the total DNA (Lo, Y.M. et al. (1998) Am. J. Hum. Genet. 62:768-775) However, recent studies suggest that fetal DNA can reach the amount of 10-20% of total DNA in the maternal circulation (Lun, F.M. et al. (2008) Clin. Chem. 54:1664-1672). In aneuploidies, one of the chromosomes is present with additional or fewer copies. For example in trisomy 21 cases, chromosome 21 is present in three copies instead of two. Therefore, the ability to distinguish normal cases from trisomy 21 cases depends on the ability to detect the extra copy of chromosome 21. However, the high levels of maternal DNA in the maternal circulation compared to the limited amount of fetal DNA further complicate quantification.

[0004] Over the last decade a large number of different methods have been applied towards the discrimination of ffDNA from circulating maternal DNA or towards ffDNA enrichment (Chan, K.C. et al. (2004) Clin. Chem. 50:88-92; Papageorgiou, E.A. et al. (2009) Am. J. Pathol. 174:1609-1618). These include DNA-based approaches, such as sequencing approaches

(Chiu, R. W. et al.(2008) Proc. Natl. Acad. Sci. USA 105:20458-20463; Fan, H.C. et al.(2008) Proc. Natl. Acad. Sci. USA 105:16266-16271) or epigenetic based approaches which focus on the investigation of the methylation status of fetal DNA either using sodium bisulfite DNA treatment (Chim, S.S. et al. (2005) Proc. Natl. Acad. Sci. USA 102:14753-14758; PCT Publication WO 2003/020974; PCT Publication WO 2005/028674), methylation-sensitive restriction enzymes (Old, R.W. et al. (2007) Reprod. Biomed. Online 15:227-235; PCT Publication WO 2005/035725) or antibodies specific to the 5-methylcytosine residues of CpG dinucleotides across the genome (Papageorgiou, E.A. et al. (2009) Am. J. Pathol. 174:1609-1618, Papageorgiou, E.A. et al. (2011) Nature Medicine 17:510-513; Tsaliki, E. et al. (2012) Prenat. Diagn. 32:996-1001; PCT Publication WO 2011/092592). Alternative approaches have targeted fetal-specific mRNA (Ng, E.K. et al. (2003) Proc. Natl. Acad. Sci. USA 100:4748-4753) or have focused on the investigation of fetal-specific proteins (Avent, N.D. et al. (2008) Semin. Fetal Neonatal Med. 13:91-98).

[0005] The implementation of next generation sequencing (NGS) technologies in the development of NIPT of aneuploidies has revolutionized the field. In 2008, two independent groups demonstrated that NIPT of trisomy 21 could be achieved using next generation massively parallel shotgun sequencing (MPSS) (Chiu, R. W. et al.(2008) Proc. Natl. Acad. Sci. USA 105:20458-20463; Fan, H.C. et al.(2008) Proc. Natl. Acad. Sci. USA 105:16266-162710). The new era of NIPT for aneuploidies has opened new possibilities for the implementation of these technologies into clinical practice. Biotechnology companies that are partly or wholly dedicated to the development of NIPT tests have initiated large scale clinical studies towards their implementation (Palomaki, G.E. et al. (2011) Genet. Med. 13:913-920; Ehrich, M. et al. (2011) Am. J. Obstet. Gynecol. 204:205e1-11; Chen, E.Z. et al. (2011) PLoS One 6:e21791; Sehnert, A.J. et al. (2011) Clin. Chem. 57:1042-1049; Palomaki, G.E. et al. (2012); Genet. Med. 14:296-305; Bianchi, D.W. et al. (2012) Obstet. Gynecol. 119:890-901; Zimmerman, B. et al. (2012) Prenat. Diagn. 32:1233-1241; Nicolaides, K.H. et al. (2013) Prenat. Diagn. 33:575-579; Sparks, A.B. et al. (2012) Prenat. Diagn. 32:3-9). Currently four companies in the United States (SEQUENOM Inc., Verinata Health, Inc., Natera and Ariosa) are offering NIPT testing using next generation sequencing approaches.

[0006] Initial NIPT approaches used massively parallel shotgun sequencing (MPSS) NGS methodologies (see e.g., US Patent No. 7,888,017; US Patent No. 8,008,018; US Patent No. 8,195,415; US Patent No. 8,296,076; US Patent No. 8,682,594; US Patent Publication 20110201507; US Patent Publication 20120270739). Thus, these approaches are whole genome-based, in which the entire maternal sample containing both maternal DNA and free fetal DNA is subjected to amplification, sequencing and analysis.

[0007] More recently, targeted-based NGS approaches for NIPT, in which only specific sequences of interest are sequenced, have been developed. For example, a SNP-based NGS approach involving targeted amplification and analysis of SNPs on chromosomes 13, 18, 21, X and Y in a single reaction has been described (Zimmerman, B. et al. (2012) Prenat. Diagn. 32:1233-1241; Nicolaides, K.H. et al. (2013) Prenat. Diagn. 33:575-579; PCT Publication WO 2011/041485; US Patent No. 8,825,412). Furthermore, an NGS-based approach has been

developed in which only specific regions of interest are sequenced wherein three probes per targeted locus are hybridized to the complementary template. Once the three probes are hybridized they are ligated and form one continuous longer probe that is then amplified and sequenced (Sparks, A.B. et al. (2012) Prenat. Diagn. 32:3-9; US Patent Publication 20120034603). The samples are analyzed using a highly multiplexed assay termed Digital ANalysis of Selected Regions (DANSR). Such targeted approaches require significantly less sequencing than the MPSS approaches, since sequencing is only performed on specific loci on the chromosome of interest rather than across the whole genome.

[0008] Additional methodologies for NGS-based approaches to NIPT are still needed, in particular approaches that can target specific sequences of interest, thereby greatly reducing the amount of sequencing needed as compared to whole genome-based approaches.

[0009] Pools of double-stranded sequences which are 100-260 bp in length and have a GC content of 19-50% are known from Parayre et al. ("Easy DNA extraction method and optimisation of PCR-Temporal Temperature Gel Electrophoresis to identify the predominant high and low GC-content bacteria from dairy products", Journal of Microbiological Methods 69 (2007) 431-441) and Li et al. ("GC Content-Based Pan-Pox Universal PCR Assays for Poxvirus Detection", Journal of Clinical Microbiology, 2010, 48(1):268-276, Epub 2009). Parayre *et al.* discloses a method to identify bacteria with low and high GC-content from dairy products, the method involving a PCR used to generate an amplicon of a DNA sequence of interest. Li *et al.* discloses PCR assays used to amplify regions of interest from the genome of poxviruses characterized by different GC content amongst different members of the chordopoxvirus family, in order to distinguish the diverse species. In contrast, the present invention relates to a kit for performing a method of testing the risk of a chromosomal abnormality in fetal DNA in a mixed sample of maternal and fetal DNA.

Summary of the Invention

[0010] The invention relates to a kit for performing a method of testing the risk of a chromosomal abnormality in fetal DNA in a mixed sample of maternal and fetal DNA, wherein the kit comprises components for:

1. (a) preparing a sequencing library from the mixed sample;
2. (b) hybridizing the sequencing library to a pool of double-stranded Target Capture Sequences (TACS), wherein the pool of TACS comprises sequences that bind to one or more chromosomes of interest comprising a chromosomal abnormality, and wherein:
 1. (i) each sequence within the pool is between 100-260 base pairs in length;
 2. (ii) each sequence within the pool binds to the chromosome(s) of interest at least 150 base pairs away, on both the 5' end and the 3' end, from regions harboring Copy Number Variations (CNVs), Segmental duplications or repetitive DNA elements; and
 3. (iii) the GC content of the TACS is between 19% and 50%;

3. (c) isolating members of the sequencing library that bind to the TACS to obtain an enriched library;
4. (d) amplifying and sequencing the enriched library; and
5. (e) performing statistical analysis on the enriched library sequences to thereby determine a risk of the chromosomal abnormality in the fetal DNA,

wherein the kit comprises a container comprising the pool of TACS, wherein the pool of TACS comprises double-stranded sequences that bind to one or more chromosomes of interest selected from the group of chromosomes 13, 18, 21, X and Y, the chromosomes comprising a chromosomal abnormality, and wherein:

1. (1) each sequence within the pool is between 100-260 base pairs in length;
2. (2) each sequence within the pool binds to the chromosome(s) of interest at least 150 base pairs away, on both the 5' end and the 3' end, from regions harboring Copy Number Variations (CNVs), Segmental duplications or repetitive DNA elements; and
3. (3) the GC content of the TACS is between 19% and 50 %.

[0011] Also disclosed but not part of the invention are methods for non-invasive prenatal testing that allow for detecting risk of chromosomal abnormalities and utilizes a targeted approach to enrich for sequences of interest prior to massive parallel sequencing and a statistical analysis approach that allows for highly accurate counting and assessment of the chromosomal constituents of maternal plasma across regions of interest. Thus, the disclosed methods reduce the amount of sequencing needed for massive parallel sequencing, allow for high throughput application, with reduced cost and a very high degree of accuracy. The disclosed methods utilize a pool of TArget Capture Sequences (TACS) to enrich for sequences of interest in a mixed sample containing both maternal and fetal DNA. In particular, the pool of TACS is designed such that the sequences within the pool have features that optimize the efficiency, specificity and accuracy of the chromosomal abnormality assessment. More specifically, the size of the TACS, the number of TACS, their placement on the chromosome(s) of interest and their GC content all have been optimized. Hybridization of the TACS to a sequencing library prepared from a mixed sample of maternal and fetal DNA (e.g., a maternal plasma sample containing ffDNA), followed by isolation of the sequences within the library that bind to the TACS allows for enrichment of only those chromosomal regions of interest, prior to massive parallel sequencing and analysis.

[0012] Also disclosed but not part of the invention is a method of testing for risk of a chromosomal abnormality in a chromosome of interest in fetal DNA in a mixed sample of maternal and fetal DNA, the method comprising:

1. (a) preparing a sequencing library from the mixed sample;
2. (b) hybridizing the sequencing library to a pool of TArget Capture Sequences (TACS), wherein the pool of TACS comprises sequences that bind to one or more chromosomes of interest and wherein:

1. (i) each sequence within the pool is between 100-260 base pairs in length;
2. (ii) each sequence within the pool binds to the chromosome(s) of interest at least 150 base pairs away, on both the 5' end and the 3' end, from regions harboring Copy Number Variations, Segmental duplications or repetitive DNA elements; and
3. (iii) the GC content of the TACS is between 19%-50%;
3. (c) isolating members of the sequencing library that bind to the TACS to obtain an enriched library;
4. (d) amplifying and sequencing the enriched library; and
5. (e) performing statistical analysis on the sequencing output of the enriched library sequences to thereby determine a risk of the chromosomal and/or other genetic abnormality in the fetal DNA.

[0013] In one embodiment, the chromosomal abnormality is an aneuploidy, such as a trisomy. The chromosome of interest can be any chromosome, although preferred chromosomes include chromosomes 13, 18, 21, X and Y. A preferred aneuploidy for detection is trisomy 21 (T21). In addition to numerical abnormalities such as aneuploidies, the invention allows for detection of other types of chromosomal abnormalities, such as structural abnormalities, including but not limited to copy number changes including, but not limited to, microdeletions and microduplications, insertions, translocations, inversions and small-size mutations including point mutations and mutational signatures

[0014] In one embodiment, the pool of TACS is fixed to a solid support. For example, the TACS can be biotinylated and bound to streptavidin-coated magnetic beads. In another embodiment, the pool of TACS may be free-moving in solution.

[0015] In one embodiment, the TACS are designed to bind to a chromosome of interest and one or more reference sequences to detect risk of a chromosomal abnormality on the chromosome of interest. Alternatively, the pool of TACS can be designed to bind to multiple chromosomes of interest such that risk of multiple chromosomal abnormalities can be detected, as well as, for example, fetal gender, all within a single analysis of the sample. For example, in one embodiment, the pool of TACS comprises of different sequences that bind to chromosomes 13, 18, 21 and X, or to chromosomes 13, 18, 21, X and Y.

[0016] In various embodiments, the GC content of the TACS is between 19% and 50%, between 19% and 49%, between 19% and 48%, between 19% and 47%, between 19% and 46%, between 19% and 45%, between 19% and 44%, between 19% and 43%, between 19% and 42%, between 19% and 41% or between 19% and 40%.

[0017] In various embodiments, each sequence within the pool of TACS is between 100-260 base pairs, 150-260 base pairs, 100-200 base pairs or 200-260 base pairs in length. In one embodiment, each sequence within the pool of TACS is 250 base pairs in length.

[0018] In various embodiments, the pool of TACS can comprise 800 or more, 1500 or more distinct sequences, 2000 or more distinct sequences, 2500 or more distinct sequences, or 3000 or more distinct sequences. In one embodiment, the pool of TACS comprises 1600 distinct sequences.

[0019] In one embodiment, sequencing of the enriched library provides a read-depth for the loci found on the chromosome of interest and read-depths for reference loci and the statistical analysis is performed by applying an algorithm that tests sequentially the read-depth of the loci found on the chromosome of interest against the read-depth of the reference loci, whereby detected differences can indicate the presence of genetic variants. The algorithm steps can include, but are not limited to: (a) removal of inadequately sequenced loci; (b) GC-content bias alleviation; and (c) ploidy status classification. In one embodiment, GC-content bias is alleviated by grouping together loci of matching GC-content.

[0020] In another embodiment, sequencing of the enriched library provides the size of fragments of cell-free genetic material captured by TACS and the statistical analysis comprises use of an algorithm that compares and contrasts the distribution of fragment-sizes from test-loci and reference-loci, whereby differences in the distribution indicate the presence of genetic variants. The algorithm steps can include, but are not limited to: (a) removal of fragment-size outliers; (b) creation of a binary distribution of fragment-sizes; and (c) testing the binary distribution of fragment-sizes originating from the region of interest against the respective distribution of the reference loci in order to classify ploidy status.

[0021] Typically, ploidy status classification is achieved by application of one or more statistical methods. For example, the statistical method can be selected from the group consisting of a t-test, a bivariate nonparametric bootstrap test, a stratified permutation test and a binomial test of proportions and/or combinations thereof. In one embodiment, all four of the aforementioned statistical methods are applied to the sample. Typically, the statistical method results in a score value for the mixed sample and risk of the chromosomal abnormality in the fetal DNA is detected when the score value for the mixed sample is above a reference threshold value. The disclosed method can further comprise estimation of the fetal DNA fraction in the mixed sample.

[0022] In another embodiment, the statistical method can be selected from the group consisting of a t-test, a bivariate nonparametric bootstrap test and a stratified permutation test. In one embodiment, all of the aforementioned statistical methods are applied to the sample. Typically, the statistical method results in a score value for the mixed sample and risk of the chromosomal abnormality in the fetal DNA is detected when the score value for the mixed sample is above a reference threshold value. The disclosed method can further comprise estimation of the fetal DNA fraction in the mixed sample.

[0023] The invention provides kits for performing the disclosed method. In one embodiment, the kit comprises a container comprising the pool of TACS and instructions for performing the method. In various other embodiments, the kit comprises additional components for carrying

out the other steps of the method.

Brief Description of the Figures

[0024]

Figure 1 is a schematic diagram of multiplexed parallel analysis of targeted genomic regions for non-invasive prenatal testing using TArget Capture Sequences (TACS).

Figure 2 is a listing of exemplary chromosomal regions for amplifying TACS that bind to chromosomes 13, 18, 21 or X.

Figure 3 is a graph of the score value assignment of 98 maternal blood samples subjected to multiplexed parallel analysis of targeted genomic regions using TACS for trisomy 21 risk detection, wherein the differences in median read-depth of the conditionally paired groups are tested for statistical significance using a t-test formula (referred to herein as statistical method 1).

Figure 4 is a graph of the score value assignment of 98 maternal blood samples subjected to multiplexed parallel analysis of targeted genomic regions using TACS for trisomy 21 risk detection, analyzed using a bivariate nonparametric bootstrap method (referred to herein as statistical method 2).

Figure 5 is a graph of the score value assignment of 98 maternal blood samples subjected to multiplexed parallel analysis of targeted genomic regions using TACS for trisomy 21 risk detection, analyzed using a stratified permutation test (referred to herein as statistical method 3).

Figure 6 is a graph of the score value assigned to 98 maternal blood samples subjected to multiplexed parallel analysis of targeted genomic regions for trisomy 21 risk detection, analyzed using a binomial test of proportions on fragment sizes (referred to herein as statistical method 4).

Figure 7 is a graph of the weighted score values of the 98 maternal blood samples resulting from analysis using statistical methods 1, 2, 3 and the weighted score method 1 as shown in Figures 3-5.

Figure 8 is a graph of the weighted score values of the 98 maternal blood samples resulting from analysis using statistical methods 1, 2, 3, 4 and the weighted score method 1 as shown in Figures 3-6.

Figure 9 is a graph of an alternative weighted approach, weighted score method 2, on the score values of the 98 maternal blood samples resulting from analysis using statistical methods 1, 2 and 3 as shown in Figures 3-5.

Figure 10 is a graph of the score value assignment of 9 synthetic samples subjected to

multiplexed parallel analysis of targeted genomic regions using TACS for 7q11.23, analyzed using an embodiment of statistical method 1 for the detection of microdeletions.

Figure 11 is a graph of the score value assignment of 9 synthetic samples subjected to multiplexed parallel analysis of targeted genomic regions using TACS for 7q11.23, analyzed using an embodiment of statistical method 1 for the detection of microduplications.

Detailed Description

[0025] The invention relates to a kit for performing a method of testing the risk of a chromosomal abnormality in fetal DNA in a mixed sample of maternal and fetal DNA, wherein the kit comprises components for:

1. (a) preparing a sequencing library from the mixed sample;
2. (b) hybridizing the sequencing library to a pool of double-stranded Target Capture Sequences (TACS), wherein the pool of TACS comprises sequences that bind to one or more chromosomes of interest comprising a chromosomal abnormality, and wherein:
 1. (i) each sequence within the pool is between 100-260 base pairs in length;
 2. (ii) each sequence within the pool binds to the chromosome(s) of interest at least 150 base pairs away, on both the 5' end and the 3' end, from regions harboring Copy Number Variations (CNVs), Segmental duplications or repetitive DNA elements; and
 3. (iii) the GC content of the TACS is between 19% and 50%;
3. (c) isolating members of the sequencing library that bind to the TACS to obtain an enriched library;
4. (d) amplifying and sequencing the enriched library; and
5. (e) performing statistical analysis on the enriched library sequences to thereby determine a risk of the chromosomal abnormality in the fetal DNA,

wherein the kit comprises a container comprising the pool of TACS, wherein the pool of TACS comprises double-stranded sequences that bind to one or more chromosomes of interest selected from the group of chromosomes 13, 18, 21, X and Y, the chromosomes comprising a chromosomal abnormality, and wherein:

1. (1) each sequence within the pool is between 100-260 base pairs in length;
2. (2) each sequence within the pool binds to the chromosome(s) of interest at least 150 base pairs away, on both the 5' end and the 3' end, from regions harboring Copy Number Variations (CNVs), Segmental duplications or repetitive DNA elements; and
3. (3) the GC content of the TACS is between 19% and 50 %.

[0026] The patent discloses a NIPT method that involves hybridization-based enrichment of

selected target regions across the human genome in a multiplexed panel assay, followed by quantification, coupled with a novel bioinformatics and mathematical analysis pipeline. In-solution hybridization enrichment has been used in the past to enrich specific regions of interest prior to sequencing (see e.g., Meyer, M and Kirchner, M. (2010) Cold Spring Harb. Protoc. 2010(6):pdbprot5448; Liao, G.J. et al. (2012) PLoS One 7:e38154; Maricic, T. et al. (2010) PLoS One 5:e14004; Tewhey, R. et al.(2009) Genome Biol. 10:R116; Tsangaras, K. et al. (2014) PLoS One 9:e109101). However, for the disclosed NIPT methods, the target sequences used to enrich for specific regions of interest relevant for detecting risk of a chromosomal abnormality have been optimized for maximum efficiency, specificity and accuracy. The human genome is full of elements that can confound and perplex any type of genetic analysis, thereby evidencing the benefit of a targeted approach for NIPT. Given this, the complexity of the human genome and the presence of these confounding elements, requires careful design of the target-capture sequences used for enrichment. As described herein, optimal Target Capture Sequences (TACS) have now been designed that allow for simpler and more robust NIPT while minimizing the risks of false positive and false negative results, that are associated with whole genome NIPT tests due to inevitable sequencing of confounding elements.

[0027] The disclosed method for testing for risk of a chromosomal abnormality in a chromosome of interest in fetal DNA in a mixed sample of maternal and fetal DNA, comprises:

1. (a) preparing a sequencing library from the mixed sample;
2. (b) hybridizing the sequencing library to a pool of Target Capture Sequences (TACS), wherein the pool of TACS comprises sequences that bind to one or more chromosomes of interest and wherein:
 1. (i) each sequence within the pool is between 100-260 base pairs in length, each sequence having a 5' end and a 3' end;
 2. (ii) each sequence within the pool binds to the chromosome(s) of interest at least 150 base pairs away, on both the 5' end and the 3' end, from regions harboring Copy Number Variations (CNVs), Segmental duplications or repetitive DNA elements; and
 3. (iii) the GC content of the TACS is between 19% and 50%;
3. (c) isolating members of the sequencing library that bind to the TACS to obtain an enriched library;
4. (d) amplifying and sequencing the enriched library; and
5. (e) performing statistical analysis on the sequencing output of the enriched library sequences to thereby determine a risk of the chromosomal abnormality in the fetal DNA.

[0028] Kits for performing the disclosed method are also encompassed.

[0029] Various aspects of this disclosure are described in further detail in the following subsections.

Target Capture Sequence Design

[0030] As used herein, the term "Target Capture Sequences" or "TACS" refers to short DNA sequences that are complementary to the region(s) of interest on a chromosome(s) of interest and which are used as "bait" to capture and enrich the region of interest from a large library of sequences, such as a whole genomic sequencing library prepared from a maternal plasma sample. A pool of TACS is used for enrichment, wherein the sequences within the pool have been optimized with regard to: (i) the length of the sequences; (ii) the distribution of the TACS across the region(s) of interest; and (iii) the GC content of the TACS. The number of sequences within the TACS pool (pool size) has also been optimized.

[0031] It has been discovered that TACS having a length of 100-260 base pairs are optimal to maximize enrichment efficiency. In various other embodiments, each sequence within the pool of TACS is between 150-260 base pairs, 100-200 base pairs, or 200-260 base pairs in length. In preferred embodiments, the length of the TACS within the pool is 250 or 260 base pairs. It will be appreciated by the ordinarily skilled artisan that a slight variation in TACS size typically can be used without altering the results (e.g., the addition or deletion of a few base pairs on either end of the TACS); accordingly, the base pair lengths given herein are to be considered "about" or "approximate", allowing for some slight variation (e.g., 1-5%) in length. Thus, for example, a length of "250 base pairs" is intended to refer to "about 250 base pairs" or "approximately 250 base pairs", such that, for example, 248 or 252 base pairs is also encompassed.

[0032] The distribution of the TACS across each region or chromosome of interest has been optimized to avoid high copy repeats, low copy repeats and copy number variants, while at the same time also being able to target informative single nucleotide polymorphisms (SNPs) in order to enable both aneuploidy, or structural copy number change detection, and fetal fraction (ff) estimation. Accordingly, each sequence within the TACS pool is designed such that the 5' end and the 3' end are each at least 150 base pairs away from regions in the genome that are known to harbour one or more of the following genomic elements: Copy Number Variations (CNVs), Segmental duplications and/or repetitive DNA elements (such as transposable elements or tandem repeat areas). In various other embodiments, each sequence within the TACS pool is designed such that the 5' end and the 3' end are each at least 200, 250, 300, 400 or 500 base pairs away from regions in the genome that are known to harbour one or more of the aforementioned elements.

[0033] The term "Copy Number Variations" is a term of art that refers to a form of structural variation in the human genome in which there can be alterations in the DNA of the genome in different individuals that can result in a fewer or greater than normal number of a section(s) of the genome in certain individuals. CNVs correspond to relatively large regions of the genome that may be deleted (e.g., a section that normally is A-B-C-D can be A-B-D) or may be duplicated (e.g., a section that normally is A-B-C-D can be A-B-C-C-D). CNVs account for

roughly 13% of the human genome, with each variation ranging in size from about 1 kilobase to several megabases in size.

[0034] The term "Segmental duplications" (also known as "low-copy repeats") is also a term of art that refers to blocks of DNA that range from about 1 to 400 kilobases in length that occur at more than one site within the genome and typically share a high level (greater than 90%) of sequence identity. Segmental duplications are reviewed in, for example, Eichler. E.E. (2001) Trends Genet. 17:661-669.

[0035] The term "repetitive DNA elements" (also known as "repeat DNA" or "repeated DNA") is also a term of art that refers to patterns of DNA that occur in multiple copies throughout the genome. The term "repetitive DNA element" encompasses terminal repeats, tandem repeats and interspersed repeats, including transposable elements. Repetitive DNA elements in NGS is discussed further in, for example, Todd, J. et al. (2012) Nature Reviews Genet. 13:36-46.

[0036] The TACS are designed with specific GC content characteristics in order to minimize data GC bias and to allow a custom and innovative data analysis pipeline. It has been determined that TACS with a GC content of 19-50% achieve optimal enrichment and perform best with cell free fetal DNA. Within the pool of TACS, different sequences can have different % GC content, although to be selected for inclusion with the pool, the % GC content of each sequence is chosen as between 19-50%. In some instances, the pool of TACS may be chosen so as to define a different % GC content range, deemed to be more suitable for the assessment of specific genetic abnormalities. Non-limiting examples of various % GC content ranges, can be between 19 and 50%, between 19% and 49%, between 19% and 48%, between 19% and 47%, between 19% and 46%, between 19% and 45%, between 19% and 44%, between 19% and 43%, between 19% and 42%, between 19% and 41% or between 19% and 40%.

As described in further detail below with respect to one embodiment of the data analysis, following amplification and sequencing of the enriched sequences, the test loci and reference loci can then be "matched" or grouped together according to their % GC content (e.g., test loci with a % GC content of 40% is matched with reference loci with a % GC content of 40%). It is appreciated that the % GC content matching procedure may allow slight variation in the allowed matched % GC range. A non-limiting instance, and with reference to the previously described example in text, a test locus with % GC content of 40% could be matched with reference loci of % GC ranging from 39-41%, thereby encompassing the test locus % GC within a suitable range.

[0037] To prepare a pool of TACS having the optimized criteria set forth above with respect to size, placement within the human genome and % GC content, both manual and computerized analysis methods known in the art can be applied to the analysis of the human reference genome. In one embodiment, a semi-automatic method is implemented where regions are firstly manually designed based on the human reference genome build 19 (hg19) ensuring that the aforementioned repetitive regions are avoided and subsequently are curated for GC-content using software that computes the % GC-content of each region based on its coordinates on the

human reference genome build 19 (hg19). In another embodiment, custom-built software is used to analyse the human reference genome in order to identify suitable TACS regions which fulfil certain criteria, such as but not limited to, %GC content, proximity to repetitive regions and/or proximity to other TACS.

[0038] The number of TACS in the pool has been carefully examined and adjusted to achieve the best balance between result robustness and assay cost/throughput. The pool typically contains at least 800 or more TACS, but can include more, such as 1500 or more TACS, 2000 or more TACS or 2500 or more TACS. It has been found that an optimal number of TACS in the pool is 1600. It will be appreciated by the ordinarily skilled artisan that a slight variation in pool size typically can be used without altering the results (e.g., the addition or removal of a small number of TACS); accordingly, the number sizes of the pool given herein are to be considered "about" or "approximate", allowing for some slight variation (e.g., 1-5%) in size. Thus, for example, a pool size of "1600 sequences" is intended to refer to "about 1600 sequences" or "approximately 1600 sequences", such that, for example, 1590 or 1610 sequences is also encompassed.

[0039] In view of the foregoing, in another aspect, the patent discloses a method for preparing a pool of TACS for use in the disclosed method for detecting risk of a chromosomal and/or other genetic abnormality, wherein the method for preparing the pool of TACS comprises: selecting regions in one or more chromosomes of interest having the criteria set forth above (e.g., at least 150 base pairs away on either end from the aforementioned repetitive sequences and a GC content of between 19% and 50%), preparing primers that amplify sequences that hybridize to the selected regions, and amplifying the sequences, wherein each sequence is 100-260 base pairs in length.

Sample Collection and Preparation

[0040] The disclosed methods are performed on a mixed sample that contains both maternal and fetal DNA. Typically the sample is a maternal plasma sample, although other tissue sources that contain both maternal and fetal DNA can be used. Maternal plasma can be obtained from a peripheral whole blood sample from a pregnant woman and the plasma can be obtained by standard methods. As little as 2-4 ml of plasma is sufficient to provide suitable DNA material for analysis according to the disclosed method. Total cell free DNA can then be extracted from the sample using standard techniques, non-limiting examples of which include a Qiasymphony protocol (Qiagen) suitable for free fetal DNA isolation or any other manual or automated extraction method suitable for cell free DNA isolation.

[0041] Following isolation, the cell free DNA of the mixed sample is used for sequencing library construction to make the sample compatible with a downstream sequencing technology, such as but not limited to Illumina Next Generation Sequencing. Typically this involves ligation of adapters onto the ends of the cell free DNA fragments, followed by amplification. Sequencing library preparation kits are commercially available. A non-limiting exemplary protocol for

sequencing library preparation is described in detail in Example 1.

Enrichment by TACS Hybridization

[0042] The region(s) of interest on the chromosome(s) of interest is enriched by hybridizing the pool of TACS to the sequencing library, followed by isolation of those sequences within the sequencing library that bind to the TACS. To facilitate isolation of the desired, enriched sequences, typically the TACS sequences are modified in such a way that sequences that hybridize to the TACS can be separated from sequences that do not hybridize to the TACS. Typically, this is achieved by fixing the TACS to a solid support. This allows for physical separation of those sequences that bind the TACS from those sequences that do not bind the TACS. For example, each sequence within the pool of TACS can be labeled with biotin and the pool can then be bound to beads coated with a biotin-binding substance, such as streptavidin or avidin. In a preferred embodiment, the TACS are labeled with biotin and bound to streptavidin-coated magnetic beads. The ordinarily skilled artisan will appreciate, however, that other affinity binding systems are known in the art and can be used instead of biotin-streptavidin/avidin. For example, an antibody-based system can be used in which the TACS are labeled with an antigen and then bound to antibody-coated beads. Moreover, the TACS can incorporate on one end a sequence tag and can be bound to a solid support via a complementary sequence on the solid support that hybridizes to the sequence tag. Furthermore in addition to magnetic beads, other types of solid supports can be used, such as polymer beads and the like.

[0043] Following enrichment of the sequence(s) of interest using the TACS, thereby forming an enriched library, the members of the enriched library are eluted from the solid support and are amplified and sequenced using standard methods known in the art. Standard Illumina Next Generation Sequencing is typically used, although other sequencing technologies can also be employed, which provides very accurate counting in addition to sequence information. To detect genetic abnormalities, such as but not limited to, aneuploidies or structural copy number changes requires very accurate counting and NGS is a type of technology that enables very accurate counting. Accordingly, for the detection of genetic abnormalities, such as but not limited to, aneuploidies or structural copy number changes, other accurate counting methods, such as digital PCR and microarrays can also be used instead of NGS. Non-limiting exemplary protocols for amplification and sequencing of the enriched library are described in detail in Example 3.

Data Analysis

[0044] The information obtained from sequencing of the enriched library is analyzed using an innovative biomathematical/biostatistical data analysis pipeline. This analysis pipeline exploits the characteristics of the TACS, and the high-efficiency of the target capture enables efficient

detection of aneuploidies or structural copy number changes, as well as other types of genetic abnormalities. Details of an exemplary analysis are described in depth in Example 4. In the analysis, first the sample's sequenced DNA fragments are aligned to the human reference genome. QC metrics are used to inspect the aligned sample's properties and decide whether the sample is suitable to undergo classification. These QC metrics can include, but are not limited to, analysis of the enrichment patterns of the loci of interest, such as for example the overall sequencing depth of the sample, the on-target sequencing output of the sample, TACS performance, GC bias expectation and fetal fraction quantification. For determining the risk of a chromosomal abnormality in the fetal DNA of the sample, an innovative algorithm is applied. The steps of the algorithm include, but are not limited to, removal of inadequately sequenced loci, read-depth and fragment-size information extraction at TACS-specific coordinates, genetic (GC-content) bias alleviation and ploidy status classification.

[0045] Ploidy status determination is achieved using one or more statistical methods, non-limiting examples of which include a t-test method, a bootstrap method, a permutation test and/or a binomial test of proportions and/or combinations thereof. It will be appreciated by the ordinarily skilled artisan that the selection and application of tests to be included in ploidy status determination is based on the number of data points available. As such, the suitability of each test is determined by various factors such as, but not limited to, the number of TACS utilized and the respective application for GC bias alleviation, if applicable. Thus, the aforementioned methods are to be taken as examples of the types of statistical analysis that may be employed and are not the only methods suitable for the determination of ploidy status. Typically, the statistical method results in a score value for the mixed sample and risk of the chromosomal abnormality in the fetal DNA is detected when the score value for the mixed sample is above a reference threshold value.

[0046] In particular, one aspect of the statistical analysis involves quantifying and alleviating GC-content bias. In addition to the challenge of detecting small signal changes in fetal DNA in the mixed sample (for example, but not limited to, additional or less genetic material from certain fetal chromosomal regions), the sequencing process itself introduces certain biases that can obscure signal detection. One such bias is the preferential sequencing/amplification of genetic regions based on their GC-content. As such, certain detection methods, such as but not limited to, read-depth based methods, need to account for such bias when examining sequencing data. Thus, the bias in the data needs to be quantified and, subsequently, suitable methods are applied to account for it such that genetic context dependencies cannot affect any statistical methods that may be used to quantify fetal genetic abnormality risk.

[0047] For example, one method of quantifying the GC-content bias is to use a locally weighted scatterplot smoothing (LOESS) technique on the sequencing data. Each targeted locus may be defined by its sequencing read-depth output and its' GC-content. A line of best fit through these two variables, for a large set of loci, provides an estimate of the expected sequencing read-depth given the GC-content. Once this GC-bias quantification step is completed, the next step is to use this information to account for possible biases in the data. One method is to normalize the read-depth of all loci by their expected read-depth (based on

each locus' GC-content). In principle, this unlinks the read-depth data from their genetic context and makes all data comparable between them. As such, data that are retrieved from different GC-content regions, such as for example, but not limited, to different chromosomes, can now be used in subsequent statistical tests for detection of any abnormalities. Thus, using the LOESS procedure, the GC bias is unlinked from the data prior to statistical testing. In one embodiment, the statistical analysis of the enriched library sequences comprises alleviating GC bias using a LOESS procedure.

[0048] In an alternative preferred embodiment, the GC-content bias is quantified and alleviated by grouping together loci of similar (matching) GC-content. Thus, conceptually this method for alleviating GC-content bias is comprised of three steps, as follows:

1. 1) identification and calculation of GC-content in the TACS;
2. 2) alleviation/accounting of GC-content bias using various matching/grouping procedures of the TACS; and
3. 3) calculation of risk of any genetic abnormalities that may be present in the fetus utilizing statistical and mathematical methods on datasets produced from step 2.

[0049] For the t-test method, the dataset is split into two groups; the test loci and the reference loci. For each group, subsets of groups are created where loci are categorized according to their GC-content as illustrated in a non-limiting example in the sample Table 1 below:

Table 1

GC	Reference loci read-depth	Test loci read-depth
40%	$x_1^{40}, x_2^{40}, \dots, x_{nx40}^{40}$	$y_1^{40}, y_2^{40}, \dots, y_{ny40}^{40}$
41%	$x_1^{41}, x_2^{41}, \dots, x_{nx41}^{41}$	$y_1^{41}, y_2^{41}, \dots, y_{ny41}^{41}$
42%	$x_1^{42}, x_2^{42}, \dots, x_{nx42}^{42}$	$y_1^{42}, y_2^{42}, \dots, y_{ny42}^{42}$
...	...	...

[0050] It is appreciated by the ordinarily skilled artisan that subgroup creation may involve encompassing a range of appropriate GC-content and/or a subset of loci that are defined by a given GC-content and/or GC-content range. Accordingly, the % GC content given in the non-limiting example of Table 1 are to be considered "about" or "approximate", allowing for some slight variation (e.g., 1-2%). Thus, for example, a % GC content of "40%" is intended to refer to "about 40%" or "approximately 40%", such that, for example, "39%-41%" GC-content loci may also be encompassed if deemed appropriate.

[0051] Hence, when referring to a particular GC-content it is understood that the reference and

test loci subgroups may comprise of any number of loci related to a particular % GC content and/or range.

[0052] Subsequently, for each GC-content subgroup, a representative read-depth is calculated. A number of methods may be utilized to choose this such as, but not limited to, the mean, median or mode of each set. Thus, two vectors of representative read-depth are created where one corresponds to the reference loci and the other to the test loci (e.g., X_m , Y_m). In one embodiment, the two vectors may be tested against each other to identify significant differences in read-depth. In another embodiment, the difference of the two vectors may be used to assess if there are significant discrepancies between the test and reference loci. The sample is attributed the score of the test.

[0053] For statistical analysis using a bootstrap approach, the dataset is split into two groups, the test loci and the reference loci. The GC-content of each locus is then calculated. Then the following procedure is performed:

A random locus is selected from the reference loci; its read-depth and GC-content are recorded. Subsequently, a random locus from the test loci is selected, with the only condition being that its' GC-content is similar to that of the reference locus. Its read-depth is recorded. It is appreciated by the ordinarily skilled artisan that GC-content similarity may encompass a range of suitable GC-content. As such, referral to a specific % GC content may be considered as "approximate" or "proximal" or "within a suitable range" (e.g 1%-2%) encompassing the specific % GC content under investigation. Thus, a reference-test locus pair of similar GC-content is created. The difference of the reference-test pair is recorded, say E_1 . The loci are then replaced to their respective groups. This process is repeated until a bootstrap sample of the same size as the number of test TACS present is created. A representative read-depth of the bootstrap sample is estimated, say E_{μ} , and recorded. A number of methods may be utilized to do so, such as but not limited to, the mean, mode or median value of the vector, and/or multiples thereof.

[0054] The process described above is repeated as many times as necessary and a distribution of E_{μ} is created. The sample is then attributed a score that corresponds to a percentile of this distribution.

[0055] For statistical analysis using a permutation test, the dataset is sorted firstly into two groups, the test-loci and the reference loci. For each group, subsets of groups are created, where loci are categorized according to their GC-content similarity (see columns 2 and 3 of the non-limiting sample Table 2 below). The number of loci present in each test subgroup is also recorded. The loci of the test group are utilized to calculate an estimate of the test-group's read-depth, say Y_{obs} . A representative number from each GC-content subgroup may be selected to do so. Any number of methods may be used to provide a read-depth estimate, such as but not limited to, the mean, median or mode of the chosen loci.

Table 2

GC	Reference loci read-depth	Test loci read-depth	test loci num	Merging of loci
40%	$x_1^{40}, x_2^{40}, \dots, x_{nx40}^{40}$	$y_1^{40}, y_2^{40}, \dots, y_{ny40}^{40}$	ny40	$x_1^{40}, \dots, x_{nx40}^{40}, y_1^{40}, \dots, y_{ny40}^{40}$
41%	$x_1^{41}, x_2^{41}, \dots, x_{nx41}^{41}$	$y_1^{41}, y_2^{41}, \dots, y_{ny41}^{41}$	ny41	$x_1^{41}, \dots, x_{nx41}^{41}, y_1^{41}, \dots, y_{ny41}^{41}$
42%	$x_1^{42}, x_2^{42}, \dots, x_{nx42}^{42}$	$y_1^{42}, y_2^{42}, \dots, y_{ny42}^{42}$	ny42	$x_1^{42}, \dots, x_{nx42}^{42}, y_1^{42}, \dots, y_{ny42}^{42}$
...	...	...	...	...

[0056] A distribution to test Yobs is then built utilizing loci irrespective of their test or reference status as follows. The test and reference loci of each GC-content subgroup (see last column of sample Table 2) are combined to allow for calculation of a new read-depth estimate. From each merged subgroup a number of loci are chosen at random, where this number is upper-bounded by the number of test-loci utilized in the original calculation of Yobs (e.g for GC content 40%, and in the context of the non-limiting sample Table 2, this number of loci may be in the range [1,ny40]). The new read-depth estimate is calculated from all the chosen loci. The procedure is iterated as many times as necessary in order to build a distribution of observed means. A sample is then attributed a score that corresponds to the position of Yobs in this distribution using a suitable transformation that accounts for the moments of the built distribution. As with the already described methods, it is appreciated that slight variation in % GC content is allowed (e.g 1%-2%), if deemed appropriate. Hence, reference to a specific GC-content could be taken as "about" or "approximate", so that for example when referring to a 40% GC-content, loci that are "approximately" or "about" 40% (e.g 39%-41%) may be utilized in the method.

[0057] For statistical analysis using a binomial test of proportions, fragment-sizes aligned to TACS-specific genomic coordinates are used. It has been shown that fragments of cell free genetic material originating from the placenta are smaller in length when compared to other cell free genetic material (Chan, K.C. (2004) Clin. Chem. 50:88-92). Hence, the statistic of interest is whether the proportion of small-size fragments aligned to a TACS-specific test-region deviates significantly from what is expected when comparing it to the respective proportion of other TACS-specific reference-regions, as this would indicate fetal genetic abnormalities.

[0058] Thus, fragment-sizes are assigned into two groups. Sizes related to the test loci are assigned to one group and fragment-sizes related to the reference loci are assigned to the other group. Subsequently, in each group, fragment sizes are distributed into two subgroups, whereby small-size fragments are assigned into one subgroup and all remaining fragments are designated to the remaining subgroup. The last step computes the proportion of small-sized fragments in each group and uses these quantities in a binomial test of proportions. The score of the test is attributed to the sample under investigation.

[0059] The final result of a sample may be given by combining one or more scores derived from the different statistical methods, non-limiting examples of which are given in Example 4.

Kits of the Invention

[0060] The invention provides kits for carrying out the disclosed methods. In one embodiment, the kit comprises a container consisting of the pool of TACS and instructions for performing the method. In one embodiment, the TACS are provided in a form that allows them to be bound to a solid support, such as biotinylated TACS. In another embodiment, the TACS are provided together with a solid support, such as biotinylated TACS provided together with streptavidin-coated magnetic beads. In various other embodiments, the kit can comprise additional components for carrying out other aspects of the method. For example, in addition to the pool of TACS, the kit can comprise one or more of the following (i) one or more components for isolating cell free DNA from a maternal plasma sample (e.g., as described in Example 1); (ii) one or more components for preparing the sequencing library (e.g., primers, adapters, linkers, restriction enzymes, ligation enzymes, polymerase enzymes and the like as described in detail in Example 1); (iii) one or more components for amplifying and/or sequencing the enriched library (e.g., as described in Example 3); and/or (iv) software for performing statistical analysis (e.g., as described in Example 4).

IV. Examples

[0061] The present invention is further illustrated by the following examples, which should not be construed as further limiting.

Example 1: Sample Collection and Library Preparation

Sample collection

[0062] Plasma samples were obtained anonymously from pregnant women after the 10th week of gestation. Protocols used for collecting samples for our study were approved by the Cyprus National Bioethics Committee, and informed consent was obtained from all participants.

Sample extraction

[0063] Cell Free DNA was extracted from 2-4ml plasma from each individual using a manual or

automated extraction method suitable for cell free DNA isolation such as for example, but not limited to, Qiasymphony protocol suitable for free fetal DNA isolation (Qiagen).

Library preparation

[0064] Extracted DNA from maternal plasma samples was used for sequencing library construction. Standard library preparation methods were used with the following modifications (Meyer, M. and Kircher, M. (2010) Cold Spring Harb. Protoc. 2010(6):pdb prot5448). A negative control extraction library was prepared separately to monitor any contamination introduced during the experiment. During this step, 5' and 3' overhangs were filled-in, by adding 12 units of T4 polymerase (NEB) while 5' phosphates were attached using 40 units of T4 polynucleotide kinase (NEB) in a 100µl reaction and subsequent incubation at 25°C for 15 minutes and then 12° C for 15 minutes. Reaction products were purified using the MinElute kit (Qiagen). Subsequently, adaptors P5 and P7 (see adaptor preparation) were ligated at 1:10 dilution to both ends of the DNA using 5 units of T4 DNA ligase (NEB) in a 40µl reaction for 20 minutes at room temperature, followed by purification using the MinElute kit (Qiagen). Nicks were removed in a fill-in reaction with 16 units of Bst polymerase (NEB) in a 40 µl reaction with subsequent incubation at 65 °C for 25 minutes and then 12 °C for 20 minutes. Products were purified using the MinElute kit (Qiagen). Library amplification was performed using a Fusion polymerase (Herculase II Fusion DNA polymerase (Agilent Technologies) or Pfusion High Fidelity Polymerase (NEB)) in 50 µl reactions and with the following cycling conditions, 95°C for 3 min; followed by 10 cycles at 95°C for 30 sec, 60°C for 30 sec, 72°C for 30 sec and finally 72°C for 3 min. The final library products were purified using the MinElute Purification Kit (Qiagen) and measured by spectrophotometry.

Adaptor preparation

[0065] Hybridization mixtures for adapter P5 and P7 were prepared (26) separately and incubated for 10sec at 95° C followed by a ramp from 95° C to 12° C at a rate of 0.1° C /sec. P5 and P7 reactions were combined to obtain a ready-to-use adapter mix (100 µM of each adapter). Hybridization mixtures were prepared as follows: P5 reaction mixture contained adaptor P5_F (500 µM) at a final concentration of 200 µM, adaptor P5+P7_R (500 µM) at a final concentration of 200 µM with 1X oligo hybridization buffer. In addition, P7 reaction mixture contained adaptor P7_F (500 µM) at a final concentration of 200 µM, adapter P5+P7_R(500 µM) at a final concentration of 200 µM with 1X oligo hybridization buffer (30). Sequences were as follows, wherein * = a phosphorothioate bond (PTO) (Integrated DNA Technologies) (Meyer, M. and Kircher, M. (2010) Cold Spring Harb. Protoc. 2010(6):pdb prot5448):

adaptor P5_F:

A*C*A*C*TCTTTCCCTACACGACGCTCTTCCG*A*T*C*T (SEQ ID NO: 1)

adaptor P7_F:

G*T*G*A*CTGGAGTTCAGACGTGTGCTCTTCCG*A*T*C*T (SEQ ID NO: 2),

adaptor_P5+P7_R:

A*G*A*T*CGGAA*G*A*G*C (SEQ ID NO: 3)

Example 2: Target Capture Sequences (TACS) Design and Preparation

[0066] Custom TACS were prepared for the detection of whole or partial chromosomal abnormalities for chromosomes 13, 18, 21, X, Y or any other chromosome, as well as other genetic abnormalities, such as but not limited to, microdeletion/microduplication syndromes, translocations, inversions, insertions, and other point or small size mutations. The genomic target-loci used for TACS design were selected based on their GC content and their distance from repetitive elements (minimum 150 bp away). TACS size can be variable. In one embodiment of the method the TACS range from 100-260 bp in size and are generated through a PCR-based approach as described below. The TACS were prepared by simplex polymerase chain reaction using standard Taq polymerase, primers designed to amplify the target-loci, and normal DNA used as template. The chromosomal regions used to design primers to amplify suitable loci on chromosomes 13, 18, 21 and X, to thereby prepare the pool of TACS for analysis of chromosomes 13, 18, 21 and X, are shown in Figure 2.

[0067] All custom TACS were generated using the following cycling conditions: 95°C for 3 min; 40 cycles at 95°C for 15 sec, 60°C for 15 sec, 72°C for 12 sec; and 72°C for 12 sec, followed by verification via agarose gel electrophoresis and purification using standard PCR clean up kits such as the Qiaquick PCR Purification Kit (Qiagen) or the NucleoSpin 96 PCR clean-up (Mackerey Nagel) or the Agencourt AMPure XP for PCR Purification (Beckman Coulter). Concentration was measured by Nanodrop (Thermo Scientific).

Example 3: TACS Hybridization and Amplification

TACS Biotinylation

[0068] TACS were prepared for hybridization, as previously described (Maricic, T. et al. (2010) PLoS One 5:e14004) with minor modifications, starting with blunt ending with the Quick Blunting Kit (NEB) and incubation at room temperature for 30 minutes. Reaction products were subsequently purified using the MinElute kit (Qiagen) and were ligated with a biotin adaptor using the Quick Ligation Kit (NEB) in a 40µl reaction at RT for 15 minutes. The reaction products were purified with the MinElute kit (Qiagen) and were denatured into single stranded DNA prior to immobilization on streptavidin coated magnetic beads (Invitrogen).

TACS Hybridization

[0069] Amplified libraries were mixed with blocking oligos (Maricic, T. *supra*) (200 μ M), 5 μ g of Cot-1 DNA (Invitrogen), 50 μ g of Salmon Sperm DNA (Invitrogen), Agilent hybridization buffer 2x, Agilent blocking agent 10X, and were heated at 95°C for 3 min to denature the DNA strands. Denaturation was followed by 30 minute incubation at 37°C to block repetitive elements and adaptor sequences. The resulting mixture was then added to the biotinylated TACS. All samples were incubated in a rotating incubator for 12- 48 hours at 66°C. After incubation, the beads were washed as described previously and DNA was eluted by heating (Maricic, T. *supra*). Eluted products were amplified using outer-bound adaptor primers. Enriched amplified products were pooled equimolarly and sequenced on an Illumina or any other suitable platform.

Example 4: Bioinformatics Sample Analysis**Human Genome Alignment**

[0070] For each sample, the bioinformatic pipeline routine described below was applied in order to align the sample's sequenced DNA fragments to the human reference genome. Targeted paired-end read fragments obtained from NGS results were processed to remove adaptor sequences and poor quality reads (Q-score<25) using the cutadapt software (Martin, M. et al. (2011) EMB.netJournal 17.1). The quality of the raw and/or processed reads as well as any descriptive statistics which aid in the assessment of quality check of the sample's sequencing output were obtained using the FastQC software (Babraham Institute (2015) *FastQC*) and/or other custom-built software. Processed reads which were at least 25 bases long were aligned to the human reference genome built hg19 (UCSC Genome Bioinformatics) using the Burrows-Wheel Alignment algorithm (Li, H. and Durbin, R. (2009) *Bioinformatics* 25:1754-1760). If relevant, duplicate reads were removed post-alignment. Where applicable, sequencing output pertaining to the same sample but processed on separate sequencing lanes, was merged to a single sequencing output file. The removal of duplicates and merging procedures were performed using the Picard tools software suite (Broad Institute (2015) *Picard*) and/or the Sambamba tools software suite (Sambamba reference, Tarasov, Artem, et al. "Sambamba: fast processing of NGS alignment formats." *Bioinformatics* 31.12 (2015): 2032-2034.).

[0071] The above software analysis resulted in a final aligned version of a sequenced sample against the human reference genome and all subsequent steps were based on this aligned version. Information in terms of Short Nucleotide Polymorphisms (SNPs) at loci of interest was obtained using bcftools from the SAMtools software suite (Li, H. et al. (2009) *Bioinformatics* 25:2078-2079). The read-depth per base, at loci of interest, was obtained using the mpileup

option of the SAMtools software suite, from here on referred to as the mpileup file. Information pertaining to the size of the aligned fragments was obtained using the view option of the SAMtools software suite, from here on referred to as the fragment-sizes file.

[0072] The mpileup file and the fragment-sizes file were processed using custom-build application programming interfaces (APIs) written in the Python and R programming languages (Python Software Foundation (2015) *Python*; The R Foundation (2015) *The R Project for Statistical Computing*). The APIs were used to determine the ploidy state of chromosomes of interest using a series of steps (collectively henceforth referred to as the "algorithm") and to also collect further descriptive statistics to be used as quality check metrics, such as but not limited to fetal fraction quantification (collectively henceforth referred to as the "QC metrics"). The APIs can also be used for the assessment of genetic abnormalities from data generated when applying the described method in cases of multiple gestation pregnancies, as well as other genetic abnormalities such as, but not limited to, microdeletions, microduplications, copy number variations, translocations, inversions, insertions, point mutations and mutational signatures.

QC Metrics

[0073] QC metrics were used to inspect an aligned sample's properties and decide whether the sample was suitable to undergo classification. These metrics were, but are not limited to:

1. (a) The enrichment of a sample. The patterns of enrichment are indicative of whether a sample has had adequate enrichment across loci of interest in a particular sequencing experiment (herein referred to as a "run"). To assess this, various metrics are assessed, non-limiting examples of which are:
 1. (i) overall sample on-target read depth,
 2. (ii) sample on-target sequencing output with respect to total mapped reads,
 3. (iii) individual TACS performance in terms of achieved read-depth,
 4. (iv) kurtosis and skewness of individual TACS enrichment and,
 5. (v) kurtosis and skewness moments that arise from all TACS.

The above checks are also taken into consideration with regards to GC-bias enrichment. Samples that fail to meet one or more of the criteria given above are flagged for further inspection, prior to classification.

2. (b) A sample's fetal fraction. Samples with an estimated fetal fraction that is below a specific threshold are not classified.

The Algorithm

[0074] The algorithm is a collection of data processing, mathematical and statistical model

routines arranged as a series of steps. The algorithm's steps aim in deciding the relative ploidy state of a chromosome of interest with respect to all other chromosomes of the sequenced sample and is used for the detection of whole or partial chromosomal abnormalities for chromosomes 13, 18, 21, X, Y or any other chromosome, as well as other genetic abnormalities such as, but not limited to, microdeletion/microduplication syndromes and other point or small size mutations. As such the algorithm can be used, but is not limited to, the detection of whole or partial chromosomal abnormalities for chromosomes 13, 18, 21, X, Y or any other chromosome, as well as other genetic abnormalities such as, but not limited to, microdeletions, microduplications, copy number variations, translocations, inversions, insertions, point mutations and other mutational signatures. The algorithm carries out, but is not limited to, two types of assessments, one pertaining to the read-depth information of each sample and the other to the distribution of fragment-sizes, across TACS-specific regions. One or more statistical tests may be associated with each type of assessment, non-limiting examples of which are given in the statistical methods described herein.

[0075] In the case of read-depth associated tests, the algorithm compares sequentially the read-depth of loci from each chromosome of interest (herein referred to as the test chromosome) against the read-depth of all other loci (herein referred to as the reference loci) to classify its ploidy state. For each sample, these steps were, but are not limited to:

1. (a) Removal of inadequately sequenced loci. The read-depth of each locus was retrieved. Loci that have not achieved a minimum number of reads, were considered as inadequately enriched and were removed prior to subsequent steps.
2. (b) Genetic (GC-content) bias alleviation. The sequencing procedure introduces discrepancies in read-depth across the loci of interest depending on their GC content. To account for such bias, a novel sequence-matching approach that increases both sensitivity and specificity to detect chromosomal aneuploidies was employed. The GC content of each locus on the test chromosome was identified and similar genetic loci were grouped together to form genetically matched groups. The procedure was repeated for the reference loci. Then, genetically matched groups from the test chromosome were conditionally paired with their genetically matched group counterparts on the reference chromosome(s). The groups may have any number of members. The conditionally matched groups were then used to assess the ploidy status of test chromosomes.
3. (c) Ploidy status determination. Ploidy status determination was achieved using a single statistical method and/or a weighted score approach on the result from the following, but not limited to, statistical methods:

Statistical Method 1: The differences in read-depth of the conditionally paired groups were tested for statistical significance using the t-test formula:

$$t = \frac{\bar{x} - \mu}{s/\sqrt{n}}$$

where t is the result of the t-test, $\bar{x}$ is the average of the differences of the conditionally paired groups, μ is the expected read-depth and is set to a value that represents insignificant read-depth differences between the two groups, s the standard deviation of

the differences of the conditionally paired groups and n the length of the vector of the conditionally paired differences. The magnitude of the t-score was then used to identify evidence, if any, against the null hypothesis of same ploidy between reference and test chromosomes. Specifically, $t \geq c_1$ (where c_1 is a predefined threshold belonging to the set of all positive numbers) shows evidence against the null. Results of the analysis of 98 maternal samples for chromosome 21 using this method are shown in Figure 3.

Statistical Method 2: Bivariate nonparametric bootstrap. The bootstrap method depends on the relationship between the random variables X (read-depth of reference loci) and Y (read-depth of test loci). Here, we treated the read depth of baits on the reference group (random variable denoted by X) as the independent covariate. The first step of the iterative procedure involved random sampling with replacement (bootstrapping) of the read-depths of loci on the reference chromosomes, i.e. $(x_1, g_1), \dots, (x_n, g_n)$, where the parameter g is known and denotes the GC-content of the chosen bait. Then, for each randomly selected reference bait (x_i, g_i) , a corresponding read depth was generated for a genetically matched locus i.e. $(y_1, g_1), \dots, (y_n, g_n)$. Thus, the bivariate data $(x_1, y_1), (x_2, y_2), \dots, (x_n, y_n)$ was arrived at, which was conditionally matched on their GC-content (parameter g_i). The differences between the read depths of the genetically matched bootstrapped values x_i and y_i were used to compute the statistic of interest in each iteration. In one embodiment this statistical measure can be, but is not limited to, the mode, mean or median of the recorded differences, and/or multiples thereof. The procedure was repeated as necessary to build up the distribution of the statistic of interest from these differences. The sample was assigned a score that corresponds to a specific percentile of the built distribution (e.g. 5th percentile). Under the null hypothesis the ploidy between chromosomes in the reference and test groups is not different. As such, samples whose score for a particular chromosome, was greater than a predefined threshold, say c_2 , were classified as statistically unlikely to have the same ploidy. Other statistical measures may be employed. Results of the analysis of 98 maternal samples for chromosome 21 using this method are shown in Figure 4.

Statistical Method 3: Stratified permutation test. The statistic of interest is the read-depth estimate of the test chromosome, denoted by $\hat{Y}_{obs}$, which is calculated using all loci of the test chromosome's genetically matched groups as follows:

$$\hat{Y}_{obs} = \frac{\sum_{j=1}^{j=T} \sum_{i=1}^{i=N_j} y_{ij}}{\sum_{j=1}^{j=T} N_j}$$

where y_{ij} is the read-depth of locus i part of the genetically matched group j (i.e loci belonging to a specific group based on their GC-content), N_j is the number of test loci part of the genetically matched group j and T the number of genetically matched groups.

[0076] Subsequently, a null distribution to test $\hat{Y}_{obs}$ was built. To do so, for each group j , the test and reference loci were combined (exchangeability under the null hypothesis), and each group j was sampled randomly up to N_j times without replacement (stratified permutation). This

created a vector of values, say y_i , and from this the vector's average value, say $\hat{y}_i$, was calculated. The procedure was repeated as necessary to build the null distribution. Finally, $\hat{Y}_{obs}$ - was studentised against the null distribution using the formula:

$$Z_{Y_{obs}} = \frac{\hat{Y}_{obs} - \hat{Y}}{\sigma_Y}$$

Where $\hat{Y}$ and σ_Y are the first and square root of the second moment of all permuted $\hat{y}_i$ statistic values. Samples whose $Z_{Y_{obs}}$ was greater than a predefined threshold, say c_3 , were statistically less likely to have the same ploidy in the reference and test groups. Results of the analysis of 98 maternal samples for chromosome 21 using this method are shown in Figure 5.

[0077] In the case of fragment-size associated tests, the algorithm computes the proportion of small-size fragments found in test-loci and compares it with the respective proportion in reference-loci as described in Statistical Method 4 below.

[0078] Statistical Method 4: Fragment Size Proportions. For each sample the number and size of fragments aligned onto the human reference genome at the corresponding TACS coordinates, is extracted. The data is subsequently filtered so as to remove fragment-sizes considered statistical outliers using the median outlier detection method. Specifically, outliers are defined as those fragments whose size is above or below the thresholds, F_{thr} , set by equation :

$$F_{thr} = F_{median} \pm (X \times IQR)$$

Where F_{median} is the median fragment-size of all fragments of a sample, X is a variable that can take values from the set of

$\mathbb{R}$

+, and IQR is the interquartile range of fragment sizes. Thereafter, a binomial test of proportions is carried out to test for supporting evidence against the null hypothesis, H_0 , where this is defined as:

H_0 : The proportion of small fragments of the test-region is not different from the proportion of small-fragments of the reference region,

In various embodiments of the invention, small fragments are defined as those fragments whose size is less than or equal to a subset of

$\mathbb{Z}$

+ that is upper-bounded by 160bp. If we define the set of all TACS as T , then the test region can be any proper subset S which defines the region under investigation, and the reference region is the relative complement of S in T . For example, in one embodiment of the invention, the set S is defined by all TACS-captured sequences of chromosome 21 and thus the reference set is defined by all TACS-captured fragments on the reference chromosomes, and/or other reference loci

The alternative hypothesis, H_1 , is defined as:

H_1 : The proportion of small fragments of the test-region is not equal to the proportion of test fragments of the reference region.

[0079] As such, and taking into account continuity correction, the following score is computed (Brown et. Al, Harrel):

$$W_{test} = (\bar{p} - p_{ref}) / \sqrt{\frac{\bar{p}(1 - \bar{p})}{N_{test}}}$$

Where

$$\bar{p} = \frac{(\bar{F} + 0.5)}{(N_{test} + 1)}$$

$$p_{ref} = \frac{(F_{ref} + 0.5)}{(N_{ref} + 1)}$$

F is the number of small-size fragments on the test-region, F_{ref} the number of small size fragments on the reference region, N_{test} the number of all fragments on the test region and N_{ref} the number of all fragments on the reference region.

[0080] For each sample, the algorithm tests sequentially the proportion of fragment sizes of regions under investigation (for example, but not limited to, chromosome 21, chromosome 18, chromosome 13) against reference regions; those not under investigation at the time of testing. For each sample a score is assigned for each test. Scores above a set-threshold, say c_4 , provide evidence against the null hypothesis. Results of the analysis of 98 maternal samples from chromosome 21 using this method are shown in Figure 6.

[0081] Weighted Score method 1: In one embodiment of the method, a weighted score was attributed to each sample s , computed as a weighted sum of all statistical methods using the formula:

$$V_s(R, F) = z_1 \max \{R_s, F_s\} + (1 - z_1) \min \{R_s, F_s\}$$

Where R_s is the run-specific corrected score arising from a weighted contribution of each read-depth related statistical method for sample s and is defined as:

$$R_s = \frac{(\sum_i w_i S_{is} - \bar{R}_r)}{\sigma_r}$$

and R_r is the run-specific median value calculated from the vector of all unadjusted read-depth related weighted scores that arise from a single sequencing run, and σ_r is a multiple of the standard deviation of R scores calculated from a reference set of 100 euploid samples. The terms $\max \{R_s, F_s\}$ and $\min \{R_s, F_s\}$ denote the maximum and minimum values of the bracketed set, respectively.

F_s is the run-specific corrected score arising from the fragment-size related statistical method and is defined as:

$$F_s = \frac{(W_{test} - \bar{R}_f)}{\sigma_f}$$

where W_{test} is as defined earlier, R_f is the run specific median calculated from the vector of all unadjusted fragment-related statistical scores that arise from a single sequencing run, and σ_f is a multiple of the standard deviation of F scores calculated from a reference set of 100 euploid samples.

[0082] A unique classification score of less than a predefined value indicates that there is no evidence from the observed data that a sample has a significant risk of aneuploidy. Results for the 98 maternal samples using the weighted score approach on a subset of the methods is shown in Figure 7, and from all methods in Figure 8.

Weighted Score method 2: In another embodiment of the method, the weighted score arising from the statistical methods described above was used to assign each sample a unique aneuploidy risk score using the formula:

$$R(t, c) = \sum_{j=0}^{j=N} w_j \frac{t_j}{c_j}$$

where R is the weighted score result, w_j the weight assigned to method j , t_j the observed score resulting from method j , and c_j the threshold of method j .

[0083] A unique classification score of less than a predefined value indicates that there is no evidence from the observed data that a sample has a significant risk of aneuploidy. Results for the 98 maternal samples using Statistical Methods 1-3 with the weighted score method 2 are shown in Figure 9.

[0084] Since all read depths from baits in the reference group were assumed to be generated from the same population, and in order to have a universal threshold, run-specific adjustments were also employed to alleviate run-specific biases.

[0085] The aforementioned method(s), are also suitable for the detection of other genetic abnormalities, such as but not limited to, subchromosomal abnormalities. A non-limiting example is the contiguous partial loss of chromosomal material leading to a state of microdeletion, or the contiguous partial gain of chromosomal material leading to a state of microduplication. A known genetic locus subject to both such abnormalities is 7q11.23. In one embodiment of statistical method 1, synthetic plasma samples of 5%, 10% and 20% fetal material were tested for increased risk of microdeletion and/or microduplication states for the genetic locus 7q11.23. Results are illustrated in Figure 10 for the cases of microdeletion and Figure 11 for the cases of microduplication.

[0086] For point mutations various binomial tests are carried out that take into consideration the fetal fraction estimate of the sample, f , the read-depth of the minor allele, r , and the total read-depth of the sequenced base, n . Two frequent, yet non-limiting examples involve assessment of the risk when the genetic abnormality is a recessive point mutation or a dominant point mutation.

[0087] In the non-limiting example of a recessive point mutation the null hypothesis tested is that both the mother and the fetus are heterozygous (minor allele frequency is 0.5) against the alternative in which the fetus is homozygous (minor allele frequency is $0.5-f/2$). A small p-value from the corresponding likelihood ratio test would indicate evidence against the null. In the non-limiting example of a dominant point mutation the null hypothesis tested is that the mother and fetus are homozygous at the given position against the alternative in which only the fetus

is heterozygous for the given position. A small p-value from the corresponding likelihood ratio test would indicate evidence against the null.

[0088] In addition to the above, fetal sex determination methods were also developed, with non-limiting examples given below. In one embodiment of the invention, fetal sex was assigned to a sample using a Poisson test using the formula:

$$Pr(r_y \leq k) = e^{-\lambda} \sum_{i=0}^{i=k} \frac{\lambda^i}{i!}$$

where

$$\lambda = \frac{fB\mu}{2}$$

and f is the fetal fraction estimate of the sample, B is the number of target sequences on chromosome Y, μ is the read-depth of the sample and k is the sum of reads obtained from all targets B . The null hypothesis of the Poisson test was that the sample is male. A value of $Pr(r_y)$ less than a threshold c_y was considered as enough evidence to reject the null hypothesis, i.e. the sample is not male. If any of the terms for computing $Pr(r_y)$ were unavailable, then the sample's sex was classified as NA (not available).

[0089] In another embodiment of the invention, fetal sex was assigned using the average read-depth of target sequences on chromosome Y. If the average read-depth of the target-sequences was over a predefined threshold, where such threshold may be defined using other sample-specific characteristics such as read-depth and fetal-fraction estimate, the fetal sex was classified as male. If the average read-depth was below such threshold then the sample was classified as female.

Fetal Fraction Estimation

[0090] Several methods have been developed to estimate fetal fraction that can be applied to singleton and/or to multiple gestation pregnancies. As such, and dependent on the type of pregnancy, the fetal fraction estimate can be obtained from either method or as a weighted estimate from a subset and/or all developed methods. Some non-limiting examples are given below.

[0091] In one embodiment, a machine learning technique has been developed based on Bayesian inference to compute the posterior distribution of fetal DNA fraction using allelic counts at heterozygous loci in maternal plasma of singleton pregnancies. Three possible informative combinations of maternal/fetal genotypes were utilized within the model to identify those fetal DNA fraction values that get most of the support from the observed data.

[0092] Let f denote the fetal DNA fraction. If the mother is heterozygous at a given genomic locus, the fetal genotype can be either heterozygous or homozygous resulting in expected minor allele frequencies at 0.5 and $0.5-f/2$, respectively. If the mother is homozygous and the

fetus is heterozygous then the expected minor allele frequency will be $f/2$. A Markov chain Monte Carlo method (a Metropolis-Hastings algorithm) (The R Foundation (2015) *The R Project for Statistical Computing*) was used with either a non-informative or an informative prior (i.e. incorporate additional information such as gestational age, maternal weight etc.) to obtain a sequence of random samples from the posterior probability distribution of fetal DNA fraction that is based on a finite mixture model.

[0093] In another embodiment, the fetal fraction estimate is computed only from the fetus-specific minor allele frequency (MAF) cluster, i.e the cluster formed when the mother is homozygous and the fetus is heterozygous for a given genomic locus. It is assumed that the mean value of the fetal fraction estimate is normally distributed as $N(2x, \sigma_x)$, where x is the mean of the fetus-specific MAF, and σ_x is the standard deviation of the fetus-specific MAF. The fetal fraction estimate is then obtained from percentiles of the computed distribution, $N(2x, \sigma_x)$.

[0094] For multiple gestation pregnancies, non-limiting examples of which include monozygotic and dizygotic twin pregnancies, triplet pregnancies and various egg and/or sperm donor cases, the fetal fraction can be estimated using information obtained from heterozygous genetic loci whose MAF value is less than a threshold, say M_{thresh} , and derived from potential fetus-specific SNPs. The ordinarily skilled artisan will appreciate that fetus specific SNPs can originate from any fetus, or from any possible combination of the fetuses or from all the fetuses of the gestation. As such, an algorithm that estimates the fetal fraction of the fetus with the smallest contribution to the total fetal content, by taking into account the combinatorial contribution of each fetus to the MAF values that define fetus-specific SNPs, and also allows for inhomogeneous contribution of fetal material to the total fetal content of plasma derived material has been developed. To this effect, a two-step approach is employed by the algorithm.

[0095] In one embodiment of the algorithm, the multiple gestation pregnancy under consideration is a dizygotic twin pregnancy. As a first step, the algorithmic implementation of the model utilizes all informative SNPs and allows for inhomogeneous fetal contribution that can be explained with a fold-difference in fetal fraction estimates of a set threshold, say cf . Specifically, if $f1$ and $f2$ represent the fetal fractions of fetus one and fetus two, and $f1 \leq f2$, then the assumption is that $f2 \leq cf \cdot f1$, with cf being a positive real number greater than or equal to 1. Under this assumption, the observed data D , defined as counts of the alternate and reference alleles at informative SNP loci, are believed to be generated from a mixture distribution of three Binomials (defined by parameters, $f1/2$, $f2/2$ and $(f1+f2)/2$), with the posterior distribution $p(f1, f2 | D)$ being proportional to the observational model which can be written as $p(f1 | f2, D) p(f2 | D)$. The posterior distribution $p(f1, f2 | D)$ is sampled with an MCMC Metropolis-Hastings algorithm using a uniform prior. The empirical quantile approach is performed on the generated data array to infer the fetal fractions.

[0096] As a second step, the algorithm runs a model-based clustering algorithm (Finite Gaussian mixture modeling fitted via EM algorithm; R-package: mclust) to identify whether there exists a separate outlier SNP cluster which is believed to be centered around $f1/2$.

Existence of such a cluster with a mean invalidating the $cf \geq f2/f1$ assumption, leads to estimation of $f1$ using only SNPs part of the identified cluster.

References

[0097]

Chris Fraley and Adrian E. Raftery (2002). Model-based Clustering, Discriminant Analysis and Density Estimation. Journal of the American Statistical Association, 97:611-631

Chris Fraley, Adrian E. Raftery, T. Brendan Murphy, and Luca Scrucca (2012). mclust Version 4 for R: Normal Mixture Modeling for Model-Based Clustering, Classification, and Density Estimation. Technical Report No. 597, Department of Statistics, University of Washington

Example 5: Results of Maternal Sample Analysis

[0098] Ninety-eight maternal samples were analyzed for trisomy chromosome 21 (T21) risk according to the methodologies described in Examples 1-4. The score values of T21 risk detection for the 98 samples using the statistical methods 1, 2, 3 and 4, are plotted in the graphs shown in Figures 3, 4, 5 and 6, respectively. Each dot represents the score value for an individual sample. The line illustrates the threshold "c" ($c=3.00$ for method 1 and $c=5.00$ for method 2, $c=4.00$ for method 3 and $c=0.91$ for method 4). Scores that exceeded the threshold line in the positive direction (i.e. score value $> c$, darker dots) were assigned a high risk of aneuploidy as opposed to scores that did not exceed the threshold (i.e. score value $< c$, lighter dots). Using all statistical methods, four samples were assigned as high risk for T21 aneuploidy. Moreover, all statistical methods indicated the same samples as being high risk even though the scores followed different distributions in the different methods.

[0099] The weighted scores that resulted from combinations of statistical methods 1, 2, 3 and 4 of the algorithm (as described in Example 4) are plotted in the graphs of Figures 7-9, wherein again each dot represents the score value for an individual sample and the line represents the threshold. The weighted score values also showed the same four samples as being high risk.

[0100] The four samples identified as being of high-risk for T21 by the above analyses were independently verified as being from pregnancies having a fetus with trisomy 21, thereby confirming the accuracy of the method for testing for risk of T21.

[0101] Furthermore, results from the analysis of microdeletion and microduplication synthetic samples are illustrated in Figures 10-11.

REFERENCES CITED IN THE DESCRIPTION

Cited references

This list of references cited by the applicant is for the reader's convenience only. It does not form part of the European patent document. Even though great care has been taken in compiling the references, errors or omissions cannot be excluded and the EPO disclaims all liability in this regard.

Patent documents cited in the description

- [US6258540B \[0002\]](#)
- [WO9107660A \[0002\]](#)
- [WO2003020974A \[0004\]](#)
- [WO2005028674A \[0004\]](#)
- [WO2005035725A \[0004\]](#)
- [WO2011092592A \[0004\]](#)
- [US7888017B \[0006\]](#)
- [US8008018B \[0006\]](#)
- [US8195415B \[0006\]](#)
- [US8296076B \[0006\]](#)
- [US8682594B \[0006\]](#)
- [US20110201507A \[0006\]](#)
- [US20120270739A \[0006\]](#)
- [WO2011041485A \[0007\]](#)
- [US8825412B \[0007\]](#)
- [US20120034603A \[0007\]](#)

Non-patent literature cited in the description

- LO, Y.M et al.Lancet, 1997, vol. 350, 485-487 [0002]
- BIANCHI, D. et al.Obstet. Gynecol., 2005, vol. 106, 841-844 [0002]
- LO, Y.M. et al.N. Engl. J. Med., 1998, vol. 339, 1734-1738 [0002]
- LO, Y.M et al.Am. J. Hum. Genet., 1998, vol. 62, 768-775 [0003]

- LUN, F.M et al.Clin. Chem., 2008, vol. 54, 1664-1672 [0003]
- CHAN, K.C. et al.Clin. Chem., 2004, vol. 50, 88-92 [0004]
- PAPAGEORGIOU, E.A. et al.Am. J. Pathol., 2009, vol. 174, 1609-1618 [0004] [0004]
- CHIU, R. W. et al.Proc. Natl. Acad. Sci. USA, 2008, vol. 105, 20458-20463 [0004] [0005]
- FAN, H.C. et al.Proc. Natl. Acad. Sci. USA, 2008, vol. 105, 16266-16271 [0004]
- CHIM, S.S. et al.Proc. Natl. Acad. Sci. USA, vol. 102, 14753-14758 [0004]
- OLD, R.W. et al.Reprod. Biomed. Online, 2007, vol. 15, 227-235 [0004]
- PAPAGEORGIOU, E.A. et al.Nature Medicine, 2011, vol. 17, 510-513 [0004]
- TSALIKI, E. et al.Prenat. Diagn., 2012, vol. 32, 996-1001 [0004]
- NG, E.K et al.Proc. Natl. Acad. Sci. USA, 2003, vol. 100, 4748-4753 [0004]
- AVENT, N.D. et al.Semin. Fetal Neonatal Med., 2008, vol. 13, 91-98 [0004]
- FAN, H.C. et al.Proc. Natl. Acad. Sci. USA, 2008, vol. 105, 16266-162710 [0005]
- PALOMAKI, G.E. et al.Genet. Med., 2011, vol. 13, 913-920 [0005]
- EHRICH, M. et al.Am. J. Obstet. Gynecol., 2011, vol. 204, 205e1-11 [0005]
- CHEN, E.Z. et al.PLoS One, 2011, vol. 6, e21791- [0005]
- SEHNERT, A.J. et al.Clin. Chem., 2011, vol. 57, 1042-1049 [0005]
- PALOMAKI, G.E. et al.Genet. Med., 2012, vol. 14, 296-305 [0005]
- BIANCHI, D.W. et al.Obstet. Gynecol., 2012, vol. 119, 890-901 [0005]
- ZIMMERMAN, B. et al.Prenat. Diag., 2012, vol. 32, 1233-1241 [0005]
- NICOLAIDES, K.H et al.Prenat. Diagn., 2013, vol. 33, 575-579 [0005] [0007]
- SPARKS, A.B. et al.Prenat. Diagn., 2012, vol. 32, 3-9 [0005] [0007]
- ZIMMERMAN, B et al.Prenat. Diag, 2012, vol. 32, 1233-1241 [0007]
- PARAYRE et al.Easy DNA extraction method and optimisation of PCR-Temporal Temperature Gel Electrophoresis to identify the predominant high and low GC-content bacteria from dairy productsJournal of Microbiological Methods, 2007, vol. 69, 431-441 [0009]
- LI et al.GC Content-Based Pan-Pox Universal PCR Assays for Poxvirus DetectionJournal of Clinical Microbiology, 2010, vol. 48, 1268-276 [0009]
- MEYER, MKIRCHNER, M.Cold Spring Harb. Protoc, 2010, vol. 6, [0026]
- LIAO, G.J. et al.PLoS One, 2012, vol. 7, e38154- [0026]
- MARICIC, T. et al.PLoS One, 2010, vol. 5, e14004- [0026] [0068]
- TEWHEY, R. et al.Genome Biol, 2009, vol. 10, R116- [0026]
- TSANGARAS, K. et al.PLoS One, 2014, vol. 9, e109101- [0026]
- EICHLER. E.E.Trends Genet., 2001, vol. 17, 661-669 [0034]
- TODD, J. et al.Nature Reviews Genet., 2012, vol. 13, 36-46 [0035]
- CHAN, K.CClin. Chem., 2004, vol. 50, 88-92 [0057]
- MEYER, M.KIRCHER, M.Cold Spring Harb. Protoc, 2010, vol. 2010, 6prot5448- [0064]
- MEYER, M.KIRCHER, M.Cold Spring Harb. Protoc., 2010, vol. 6, prot5448- [0065]
- MARTIN, M. et al.EMB.netJournal, 2011, vol. 17, 1- [0070]
- LI, H.DURBIN, R.Bioinformatics, 2009, vol. 25, 1754-1760 [0070]
- TARASOV, ARTEM et al.Sambamba: fast processing of NGS alignment formatsBioinformatics, 2015, vol. 31, 122032-2034 [0070]
- LI, H et al.Bioinformatics, 2009, vol. 25, 2078-2079 [0071]

- **CHRIS FRALEYADRIAN E. RAFTERY**Model-based Clustering, Discriminant Analysis and Density EstimationJournal of the American Statistical Association, 2002, vol. 97, 611-631 [0097]
- **CHRIS FRALEYADRIAN E. RAFTERYT. BRENDAN MURPHYLUCA SCRULLA**mclust Version 4 for R: Normal Mixture Modeling for Model-Based ClusteringClassification, and Density Estimation, 2012, [0097]

P A T E N T K R A V

1. Kit til at fremføre en fremgangsmåde til at teste risikoen for en kromosomal abnormitet i foster DNA i en blandet prøve af moderens og fosterets DNA, hvor kittet omfatter bestanddelen til:

- 5 (a) at forberede et sekventeringsbibliotek fra den blandede prøve;
- (b) at hybridisere af sekventeringsbiblioteket til en pulje af dobbeltstrengede TArget Capture Sequences (TACS), hvor puljen af TACS omfatter sekvenser, der binder til et eller flere kromosomer af interesse, der omfatter en kromosomal abnormitet, og hvor:
- (i) hver sekvens i puljen har en længde på mellem 100-260 basepar;
- 10 (ii) hver sekvens i puljen binder til kromosomet(erne) af interesse på mindst 150 basepar væk, på både 5' enden og 3' enden, fra regionerne der indeholder Copy Number Variations (CNVs), Segmental duplikationer eller gentagende DNA elementer; og
- (iii) GC indholdet af TACS er mellem 19% og 50%;
- (c) at isolere medlemmer af sekventeringsbiblioteket som binder til TACS for at op-
- 15 nå et beriget bibliotek;
- (d) at opformere og sekventere det beriget bibliotek; og
- (e) at udføre statistik analyse på det beriget bibliotek sekvenser til derved at bestemme en risiko for en kromosomal abnormitet i foster DNA,
- hvor kittet omfatter en beholder der omfatter puljen af TACS, hvor puljen af TACS
- 20 omfatter dobbelt-strengt sekvenser, som binder til én eller flere kromosomer af interesse valgt fra gruppen af kromosomerne 13, 18, 21, X og Y, kromosomerne der omfatter en kromosomal abnormitet, og hvor:
- (1) hver sekvens i puljen er mellem 100-260 basepar lang;
- (2) hver sekvens i puljen binder til kromosomet(erne) af interesse på mindst 150
- 25 basepar væk, på både 5' enden og 3' enden, fra regionerne der indeholder Copy Number Variants (CNVs), Segmental duplikationer eller gentagende DNA elementer; og
- (3) GC indholdet af TACS er mellem 19% og 50%.

DRAWINGS

Drawing




FIGURE 1

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr13	91049752	91050001	chr13	21245232	21245481	chr13	28077809	28078058
chr13	92177139	92177388	chr13	21251423	21251672	chr13	28212786	28213035
chr13	93071049	93071298	chr13	21276601	21276850	chr13	28302961	28303210
chr13	94326973	94327222	chr13	21296989	21297237	chr13	28325874	28326123
chr13	94415507	94415756	chr13	21321012	21321261	chr13	28395893	28396142
chr13	95600034	95600283	chr13	21330138	21330386	chr13	28549591	28549840
chr13	112910529	112910778	chr13	21348854	21349102	chr13	28592602	28592851
chr13	53328405	53328654	chr13	21351692	21351941	chr13	28608293	28608542
chr13	20533985	20534234	chr13	21356527	21356775	chr13	28617043	28617292
chr13	20535728	20535977	chr13	21361621	21361869	chr13	28622471	28622720
chr13	20544227	20544476	chr13	21396343	21396591	chr13	28623594	28623843
chr13	20567359	20567607	chr13	21399744	21399992	chr13	28941450	28941699
chr13	20577159	20577407	chr13	21417352	21417601	chr13	28958898	28959147
chr13	20598772	20599021	chr13	21422981	21423229	chr13	28960375	28960624
chr13	20600556	20600804	chr13	21607075	21607324	chr13	29003985	29004234
chr13	20610779	20611027	chr13	21631007	21631256	chr13	29052686	29052935
chr13	20621336	20621585	chr13	22114658	22114907	chr13	52293092	52293341
chr13	20624392	20624641	chr13	22130956	22131205	chr13	52302656	52302905
chr13	20625461	20625710	chr13	22189915	22190164	chr13	53286404	53286653
chr13	20632609	20632857	chr13	22270443	22270692	chr13	53306303	53306552
chr13	20635196	20635444	chr13	22396225	22396474	chr13	53497224	53497473
chr13	21127405	21127653	chr13	22609278	22609527	chr13	53548644	53548893
chr13	21163998	21164247	chr13	22691882	22692131	chr13	53652715	53652964
chr13	21170950	21171199	chr13	22765144	22765393	chr13	53861872	53862121
chr13	21172515	21172764	chr13	27175379	27175628	chr13	53897429	53897678
chr13	21173774	21174022	chr13	27405837	27406086	chr13	53938205	53938454
chr13	21182580	21182828	chr13	27646720	27646969	chr13	54626806	54627055
chr13	21203385	21203634	chr13	27674747	27674996	chr13	54740174	54740423
chr13	21217173	21217422	chr13	27969360	27969609	chr13	54753384	54754083
chr13	21228403	21228652	chr13	27978231	27978480	chr13	54768238	54768487
chr13	21237418	21237666	chr13	28003481	28003730	chr13	54773333	54773582
chr13	21239719	21239968	chr13	28073180	28073429	chr13	54910969	54911218

FIGURE 2 (Page 1 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr13	54917347	54917596	chr13	111956688	111956937	chr13	60654975	60655224
chr13	54926251	54926500	chr13	112517966	112518215	chr13	61151386	61151635
chr13	55009732	55009981	chr13	112531519	112531768	chr13	62064603	62064852
chr13	59785342	59785591	chr13	112986709	112986958	chr13	71479359	71479608
chr13	60086302	60086551	chr13	113000054	113000303	chr13	39233534	39233783
chr13	60668020	60668269	chr13	21391429	21391678	chr13	39281657	39281916
chr13	60707747	60707996	chr13	21394616	21394865	chr13	39343524	39343773
chr13	60975949	60976198	chr13	21404268	21404516	chr13	39444897	39445146
chr13	61177164	61177413	chr13	21618233	21618482	chr13	39469458	39469707
chr13	61495655	61495904	chr13	21619521	21619769	chr13	39486237	39486486
chr13	67720914	67721163	chr13	21629602	21629851	chr13	39541262	39541511
chr13	67799557	67799806	chr13	21654127	21654376	chr13	39597307	39597556
chr13	70101783	70102032	chr13	22215795	22216044	chr13	39662550	39662799
chr13	71512697	71512946	chr13	22227261	22227510	chr13	39681445	39681694
chr13	71821654	71821903	chr13	22243011	22243260	chr13	39707053	39707302
chr13	72405844	72406093	chr13	22255153	22255402	chr13	39791099	39791348
chr13	74087455	74087704	chr13	22639387	22639636	chr13	39882804	39883053
chr13	74192477	74192726	chr13	27183101	27183350	chr13	39887397	39887646
chr13	74260893	74261142	chr13	27387690	27387939	chr13	40023278	40023527
chr13	78218064	78218313	chr13	30466189	30466438	chr13	40038277	40038526
chr13	88033439	88033688	chr13	30977198	30977447	chr13	40071926	40072175
chr13	90954378	90954627	chr13	31061235	31061484	chr13	40138914	40139163
chr13	91290109	91290358	chr13	31205037	31205286	chr13	40171414	40171663
chr13	91874767	91875016	chr13	32625308	32625557	chr13	40405751	40406000
chr13	92125914	92126163	chr13	34428102	34428351	chr13	40483291	40483540
chr13	92308079	92308328	chr13	36054899	36055148	chr13	40562330	40562579
chr13	92971372	92971621	chr13	36105686	36105935	chr13	40590631	40590880
chr13	94063232	94063481	chr13	40676702	40676951	chr13	40633944	40634193
chr13	94453836	94454085	chr13	40745383	40745632	chr13	40644042	40644291
chr13	94713643	94713892	chr13	49281723	49281972	chr13	40661582	40661831
chr13	95766126	95766375	chr13	53624527	53624776	chr13	40667939	40668188
chr13	95814599	95814848	chr13	60069277	60069526	chr13	40684306	40684555

FIGURE 2 (Page 2 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr13	40747946	40748195	chr13	47314093	47314342	chr13	30145894	30146143
chr13	40767281	40767530	chr13	47623498	47623747	chr13	30187253	30187502
chr13	41165017	41165266	chr13	48562058	48562307	chr13	30421518	30421767
chr13	41176869	41177118	chr13	48698166	48698415	chr13	30531947	30532196
chr13	41573543	41573792	chr13	48762184	48762433	chr13	30778718	30778967
chr13	41581284	41581533	chr13	48806515	48806764	chr13	30905760	30906009
chr13	41588844	41589093	chr13	48811339	48811588	chr13	31011560	31011809
chr13	41628316	41628565	chr13	49119202	49119451	chr13	31089936	31090185
chr13	41633944	41634193	chr13	49479756	49480005	chr13	31227324	31227573
chr13	41827812	41828061	chr13	49505070	49505319	chr13	31430437	31430686
chr13	42280338	42280587	chr13	49530727	49530976	chr13	31575060	31575309
chr13	42349979	42350228	chr13	50387096	50387317	chr13	31641685	31641934
chr13	42589616	42589865	chr13	50960526	50960775	chr13	31718167	31718416
chr13	42773088	42773337	chr13	51783208	51783457	chr13	31723797	31724044
chr13	42811933	42812182	chr13	51816719	51816968	chr13	32249897	32250146
chr13	44211957	44212206	chr13	51825241	51825490	chr13	32296888	32297137
chr13	44239246	44239495	chr13	52025153	52025402	chr13	32311823	32312072
chr13	44247433	44247682	chr13	52249504	52249753	chr13	32316395	32316644
chr13	44359516	44359765	chr13	29111011	29111260	chr13	32321163	32321412
chr13	44388992	44389241	chr13	29192994	29193243	chr13	32348262	32348511
chr13	44493443	44493692	chr13	29233772	29234021	chr13	32565108	32565357
chr13	44557886	44558135	chr13	29310139	29310388	chr13	32578964	32579213
chr13	44959162	44959411	chr13	29563140	29563389	chr13	32584041	32584290
chr13	45043950	45044199	chr13	29580073	29580322	chr13	32596887	32597136
chr13	45115235	45115484	chr13	29929175	29929424	chr13	32651705	32651954
chr13	45147386	45147635	chr13	29934475	29934724	chr13	32760468	32760717
chr13	47067063	47067312	chr13	29939640	29939889	chr13	32829312	32829561
chr13	47120595	47120844	chr13	29957853	29958102	chr13	32852622	32852871
chr13	47154937	47155186	chr13	29995732	29995981	chr13	32872003	32872252
chr13	47258834	47259083	chr13	30003375	30003624	chr13	33535420	33535669
chr13	47261049	47262098	chr13	30014954	30015203	chr13	33634729	33634978
chr13	47286641	47286890	chr13	30114022	30114271	chr13	33660366	33660615

FIGURE 2 (Page 3 of 31)

chromosome	start	stop
chr13	33742309	33742558
chr13	33760763	33761012
chr13	33787480	33787729
chr13	33803219	33803468
chr13	33925858	33926107
chr13	34286418	34286667
chr13	34301780	34302029
chr13	34380783	34381032
chr13	35504924	35505173
chr13	35532835	35533084
chr13	35538772	35539021
chr13	36233165	36233414
chr13	36343466	36343715
chr13	36348003	36348252
chr13	36380590	36380839
chr13	36386774	36387023
chr13	36410707	36410956
chr13	36427833	36428082
chr13	364445018	36445267
chr13	36451093	36451342
chr13	36568336	36568585
chr13	36576882	36577131
chr13	36663599	36663848
chr13	37595672	37595921
chr13	38283224	38283473
chr13	38320894	38321143
chr13	38361974	38362223
chr13	38380696	38380945
chr13	38789453	38789702
chr13	39009512	39009761
chr13	39096011	39096260
chr13	39202882	39203131

chromosome	start	stop
chr13	39264689	39264941
chr13	39265199	39265451
chr13	39265501	39265748
chr13	39266515	39266762
chr13	39266915	39267157
chr13	39322100	39322347
chr13	39338010	39338257

chromosome	start	stop
------------	-------	------

FIGURE 2 (Page 4 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chrX	44004660	44004909	chrX	96475190	96475439	chrX	132160888	132161137
chrX	44005572	44005821	chrX	96522258	96522507	chrX	133560531	133560780
chrX	80622102	80622351	chrX	96605785	96606034	chrX	133847304	133847553
chrX	81576088	81576337	chrX	98226881	98227130	chrX	133886077	133886326
chrX	83458870	83459119	chrX	98227467	98227716	chrX	133886778	133887027
chrX	83461989	83462238	chrX	98314697	98314946	chrX	135279738	135279987
chrX	83972195	83972444	chrX	98331718	98331967	chrX	135581031	135581280
chrX	83973091	83973340	chrX	102986351	102986600	chrX	136941543	136941792
chrX	84308155	84308404	chrX	112117885	112118134	chrX	138913313	138913562
chrX	84412933	84413182	chrX	116461543	116461792	chrX	138913920	138914169
chrX	84530826	84531075	chrX	116746667	116746916	chrX	139857034	139857283
chrX	84531356	84531605	chrX	117121287	117121536	chrX	139857318	139857567
chrX	84561384	84561633	chrX	117494627	117494876	chrX	141400062	141400311
chrX	85583338	85583587	chrX	118400460	118400709	chrX	141408144	141408393
chrX	85584098	85584347	chrX	119760145	119760394	chrX	141408548	141408797
chrX	85643367	85643616	chrX	120436725	120436974	chrX	142517923	142518172
chrX	85643847	85644096	chrX	120437307	120437556	chrX	143466384	143466633
chrX	85644765	85645014	chrX	121854317	121854566	chrX	147766073	147766322
chrX	85682212	85682461	chrX	121855104	121855353	chrX	147869477	147869726
chrX	85715710	85715959	chrX	122863038	122863287	chrX	147880766	147881015
chrX	85764768	85765017	chrX	123203756	123204005	chrX	147881811	147882060
chrX	85988749	85988998	chrX	123989030	123989279	chrX	151796091	151796340
chrX	85989214	85989463	chrX	124394967	124395216	chrX	104102171	104102420
chrX	85997574	85997823	chrX	125455385	125455634	chrX	104124449	104124698
chrX	87225841	87226090	chrX	125500101	125500350	chrX	104234605	104234854
chrX	91517109	91517358	chrX	125500582	125500831	chrX	104235195	104235444
chrX	91517578	91517827	chrX	127973758	127974007	chrX	104291380	104291629
chrX	94601042	94601291	chrX	128274367	128274616	chrX	104291786	104292035
chrX	95282210	95282459	chrX	128274689	128274938	chrX	104647680	104647929
chrX	95700693	95700942	chrX	128275238	128275487	chrX	104648691	104648940
chrX	96016317	96016566	chrX	128694864	128695113	chrX	106886942	106887191
chrX	96474722	96474971	chrX	131201343	131201592	chrX	3343487	3343736

FIGURE 2 (Page 5 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chrX	4013011	4013260	chrX	10710087	10710336	chrX	18600682	18600931
chrX	5008870	5009119	chrX	11139863	11140112	chrX	19376384	19376633
chrX	5687782	5688031	chrX	11658912	11659161	chrX	20026004	20026253
chrX	5717003	5717252	chrX	11681358	11681607	chrX	20627108	20627357
chrX	5808998	5809247	chrX	11686196	11686445	chrX	21537202	21537451
chrX	5814187	5814436	chrX	11686913	11687162	chrX	21629967	21630216
chrX	5836574	5836823	chrX	11687182	11687431	chrX	21887062	21887311
chrX	5849776	5850025	chrX	11917847	11918096	chrX	22202443	22202692
chrX	6895474	6895723	chrX	12016474	12016723	chrX	23018625	23018874
chrX	6965226	6965475	chrX	12158815	12159064	chrX	23019006	23019255
chrX	6965573	6965822	chrX	12485152	12485401	chrX	24521525	24521774
chrX	7152743	7152992	chrX	12608017	12608266	chrX	24522047	24522296
chrX	9529739	9529988	chrX	12724960	12725209	chrX	25401534	25401783
chrX	9605323	9605572	chrX	12841791	12842040	chrX	25403300	25403549
chrX	9615158	9615407	chrX	12904556	12904807	chrX	28689766	28690015
chrX	9766329	9766578	chrX	12937788	12938037	chrX	28693600	28693849
chrX	9793372	9793621	chrX	12938100	12938349	chrX	28717026	28717275
chrX	9849298	9849547	chrX	12992266	12992515	chrX	28720131	28720380
chrX	9850182	9850431	chrX	12995926	12996175	chrX	28808128	28808377
chrX	10002581	10002830	chrX	12997273	12997522	chrX	28824813	28825062
chrX	10301512	10301761	chrX	13293562	13293811	chrX	28825844	28826093
chrX	10354762	10355011	chrX	13338207	13338456	chrX	29487025	29487274
chrX	10375361	10375610	chrX	13340422	13340671	chrX	30847194	30847443
chrX	10384260	10384509	chrX	13351532	13351781	chrX	32660941	32661190
chrX	10398645	10398894	chrX	13624857	13625106	chrX	32949975	32950224
chrX	10572719	10572968	chrX	15525635	15525884	chrX	35971182	35971431
chrX	10579800	10580049	chrX	15543469	15543718	chrX	35971544	35971793
chrX	10580482	10580731	chrX	15543853	15544102	chrX	43599991	43600240
chrX	10627502	10627751	chrX	15844810	15845059	chrX	43600262	43600511
chrX	10645064	10645313	chrX	16860721	16860970	chrX	43807492	43807741
chrX	10669664	10669913	chrX	16861729	16861978	chrX	63137429	63137678
chrX	10709674	10709923	chrX	18307677	18307926	chrX	63137928	63138177

FIGURE 2 (Page 6 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chrX	74269728	74269977	chrX	78401753	78402002	chrX	136320360	136320609
chrX	74270926	74271175	chrX	78403248	78403497	chrX	139380500	139380749
chrX	74282022	74282271	chrX	78405536	78405785	chrX	139381074	139381323
chrX	74329664	74329913	chrX	78411791	78412040	chrX	144897923	144898172
chrX	74344769	74345018	chrX	78427326	78427575	chrX	144908837	144909086
chrX	74741441	74741690	chrX	78430991	78431240	chrX	144958353	144958602
chrX	76192017	76192266	chrX	78431372	78431621	chrX	144959377	144959626
chrX	76194203	76194452	chrX	78437851	78438100	chrX	145799354	145799603
chrX	76221949	76222198	chrX	78526539	78526788	chrX	145922717	145922966
chrX	76226569	76226818	chrX	78531686	78531935	chrX	146994777	146994726
chrX	76234682	76234931	chrX	78996675	78996924	chrX	146995311	146995560
chrX	76308591	76308840	chrX	107682036	107682285	chrX	147549777	147541226
chrX	76647025	76647274	chrX	107683930	107684179	chrX	149764787	149765036
chrX	76711476	76711725	chrX	110366745	110366994	chrX	150808513	150808762
chrX	76759387	76759636	chrX	110367927	110368176	chrX	150809248	150809497
chrX	76778190	76778439	chrX	110368330	110368579	chrX	150942801	150943050
chrX	76778847	76779096	chrX	110443307	110443556	chrX	79712826	79713075
chrX	76875971	76876220	chrX	110466084	110466333	chrX	80374192	80374441
chrX	76923678	76923927	chrX	110537791	110538040	chrX	80491243	80491492
chrX	76937441	76937590	chrX	110538060	110538309	chrX	80514786	80515035
chrX	76938848	76939097	chrX	110538444	110538693	chrX	80553311	80553560
chrX	76939481	76939730	chrX	110575740	110575989	chrX	104028154	104028403
chrX	76949897	76950146	chrX	110576099	110576348	chrX	104039199	104039448
chrX	76951314	76951563	chrX	110950649	110950898	chrX	104101675	104101924
chrX	77040328	77040577	chrX	110971144	110971393			
chrX	77168344	77168593	chrX	121597454	121597703			
chrX	77634270	77634519	chrX	135558463	135558712			
chrX	78176076	78176325	chrX	135570815	135571064			
chrX	78176555	78176804	chrX	135760826	135761075			
chrX	78178901	78179150	chrX	135831981	135832230			
chrX	78208245	78208494	chrX	136044394	136044643			
chrX	78290291	78290540	chrX	136318161	136318410			

FIGURE 2 (Page 7 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr21	17021878	17022127	chr21	38751525	38751774	chr21	35093305	35093554
chr21	17046718	17046967	chr21	39173262	39173511	chr21	35094930	35095179
chr21	17088576	17088825	chr21	39285267	39285516	chr21	35099418	35099667
chr21	17134784	17135033	chr21	39330340	39330589	chr21	35138457	35138706
chr21	17161709	17161958	chr21	39442392	39442641	chr21	35152877	35153126
chr21	17183683	17183932	chr21	39499219	39499468	chr21	35157173	35157422
chr21	17301953	17302232	chr21	39541298	39541547	chr21	35161795	35162044
chr21	17302814	17303063	chr21	39563720	39563969	chr21	35296082	35296331
chr21	17336649	17336898	chr21	39584873	39585122	chr21	35450331	35450580
chr21	17346228	17346477	chr21	46530377	46530731	chr21	35473022	35473271
chr21	17348050	17348299	chr21	46534043	46534292	chr21	35473971	35474220
chr21	17351309	17351558	chr21	46535763	46536012	chr21	35498164	35498413
chr21	17362239	17362488	chr21	46537454	46537703	chr21	35540943	35541192
chr21	17367713	17367962	chr21	46542089	46542338	chr21	35702965	35703214
chr21	17506670	17506919	chr21	46556856	46557105	chr21	28574819	28575068
chr21	17548792	17549041	chr21	46564465	46564714	chr21	28586377	28586626
chr21	17591325	17591574	chr21	46567255	46567504	chr21	28597009	28597258
chr21	17695933	17696182	chr21	46572969	46573218	chr21	28616121	28616370
chr21	19087043	19087292	chr21	46605329	46605578	chr21	28631408	28631657
chr21	19178365	19178614	chr21	32134809	32135058	chr21	28891355	28891604
chr21	22579319	22579568	chr21	32148136	32148385	chr21	28892223	28892472
chr21	22605212	22605461	chr21	32165994	32166243	chr21	28919446	28919695
chr21	22964667	22964916	chr21	32174173	32174422	chr21	29012106	29012355
chr21	22984832	22985081	chr21	32179882	32180131	chr21	29021007	29021256
chr21	22985386	22985635	chr21	32181384	32181633	chr21	29753618	29753867
chr21	36076596	36076845	chr21	32184463	32184712	chr21	29796479	29796728
chr21	36169363	36169617	chr21	32201994	32202243	chr21	29812165	29812414
chr21	36789088	36789337	chr21	32208609	32208858	chr21	29813966	29814215
chr21	37939587	37939836	chr21	32319600	32319849	chr21	29814951	29815200
chr21	37959766	37960015	chr21	35034581	35034830	chr21	29820028	29820277
chr21	37983316	37983565	chr21	35060177	35060426	chr21	29825992	29826241
chr21	38559539	38559788	chr21	35070605	35070854	chr21	29827616	29827865

FIGURE 2 (Page 8 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr21	29841662	29841911	chr21	31733443	31733692	chr21	27623692	27623941
chr21	29884231	29884480	chr21	31768996	31769245	chr21	27763679	27763928
chr21	29898255	29898504	chr21	31886189	31886438	chr21	27894000	27894249
chr21	30368230	30368479	chr21	31886808	31887057	chr21	27926498	27926747
chr21	30413906	30414155	chr21	32119833	32120082	chr21	28077177	28077426
chr21	30415555	30415804	chr21	32125263	32125512	chr21	28122462	28122711
chr21	30428953	30429202	chr21	25795207	25795456	chr21	28161337	28161586
chr21	30477650	30477899	chr21	25802273	25802522	chr21	28180925	28181174
chr21	30528161	30528410	chr21	26355452	26355702	chr21	28275754	28276003
chr21	30550132	30550381	chr21	26365770	26366019	chr21	28291449	28291698
chr21	30580738	30580987	chr21	26373440	26373689	chr21	28294820	28295069
chr21	30600446	30600695	chr21	26374774	26375023	chr21	28295351	28295600
chr21	30628668	30628917	chr21	26390279	26390528	chr21	28295857	28296106
chr21	30660609	30660858	chr21	26390951	26391200	chr21	28311866	28312115
chr21	30680391	30680640	chr21	26402802	26403051	chr21	28334784	28335033
chr21	30720214	30720463	chr21	26496769	26497018	chr21	28354814	28355063
chr21	30782462	30782711	chr21	26515400	26515649	chr21	28415341	28415590
chr21	30810199	30810448	chr21	26571662	26571911	chr21	28433186	28433435
chr21	30849402	30849651	chr21	26589786	26590035	chr21	28450226	28450475
chr21	30868525	30868774	chr21	26604824	26605073	chr21	28474083	28474332
chr21	30988776	30989025	chr21	26933670	26933919	chr21	28491062	28491311
chr21	31012406	31012655	chr21	26950703	26950952	chr21	28509439	28509688
chr21	31032720	31032969	chr21	27045134	27045383	chr21	28534176	28534425
chr21	31085356	31085605	chr21	27398390	27398639	chr21	40042441	40042690
chr21	31494755	31495004	chr21	27441346	27441595	chr21	40427338	40427587
chr21	31560667	31560916	chr21	27465490	27465739	chr21	40509976	40510225
chr21	31605973	31606222	chr21	27485182	27485431	chr21	40531291	40531540
chr21	31621826	31622075	chr21	27486006	27486255	chr21	40578323	40578572
chr21	31653822	31654071	chr21	27502429	27502678	chr21	40601158	40601407
chr21	31670830	31671079	chr21	27520192	27520441	chr21	40630577	40630826
chr21	31691347	31691596	chr21	27539803	27540052	chr21	41280723	41280972
chr21	31708936	31709185	chr21	27623130	27623379	chr21	41607273	41607522

FIGURE 2 (Page 9 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr21	41629820	41630069	chr21	15843963	15844212	chr21	16265484	16265733
chr21	41703462	41703711	chr21	15858450	15858699	chr21	16273734	16273983
chr21	41729089	41729338	chr21	15859015	15859264	chr21	16275865	16276114
chr21	41761636	41761885	chr21	16002625	16002874	chr21	16290234	16290483
chr21	41839021	41839270	chr21	16030414	16030663	chr21	16295877	16296126
chr21	41908902	41909151	chr21	16063119	16063368	chr21	16296687	16296936
chr21	41933321	41933570	chr21	16068358	16068607	chr21	16298842	16299091
chr21	41981868	41982117	chr21	16071495	16071744	chr21	16316213	16316462
chr21	42029142	42029391	chr21	16081282	16081531	chr21	16319836	16320085
chr21	42030642	42030891	chr21	16115078	16115327			
chr21	42080190	42080439	chr21	16122085	16122334			
chr21	42251219	42251468	chr21	16130688	16130937			
chr21	42559559	42559808	chr21	16133865	16134114			
chr21	42601828	42602077	chr21	16139969	16140218			
chr21	42851340	42851589	chr21	16141592	16141841			
chr21	42872713	42872962	chr21	16160517	16160766			
chr21	43032728	43032977	chr21	16165861	16166110			
chr21	43406661	43406910	chr21	16189746	16189995			
chr21	46521300	46521549	chr21	16190806	16191055			
chr21	163233548	16323797	chr21	16192618	16192867			
chr21	16330249	16330498	chr21	16195809	16196058			
chr21	16332622	16332871	chr21	16198337	16198586			
chr21	16333315	16333564	chr21	16203731	16203980			
chr21	16335355	16335604	chr21	16207035	16207284			
chr21	16344361	16344610	chr21	16209341	16209590			
chr21	16349606	16349855	chr21	16211270	16211519			
chr21	16358650	16358899	chr21	16219907	16220156			
chr21	15713251	15713500	chr21	16230785	16231034			
chr21	15714357	15714606	chr21	16236026	16236275			
chr21	15725977	15726226	chr21	16241139	16241388			
chr21	15731366	15731615	chr21	16242812	16243061			
chr21	15816232	15816481	chr21	16246944	16247193			

FIGURE 2 (Page 10 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr18	64634331	64634580	chr18	73871422	73871671	chr18	50308727	50308976
chr18	64645925	64646174	chr18	73876103	73876352	chr18	50309031	50309280
chr18	64975938	64976187	chr18	73879805	73880054	chr18	50338747	50338996
chr18	64996573	64996822	chr18	73884248	73884497	chr18	50423284	50423533
chr18	65066425	65066674	chr18	73931865	73932114	chr18	50424422	50424671
chr18	65073729	65073978	chr18	73938014	73938263	chr18	50424789	50425038
chr18	65074170	65074419	chr18	73943212	73943461	chr18	50448837	50449086
chr18	66343957	66344206	chr18	73947628	73947877	chr18	50449492	50449741
chr18	66344406	66344655	chr18	73964168	73964417	chr18	50514079	50514328
chr18	66585247	66585496	chr18	73966321	73966570	chr18	50519889	50520138
chr18	66596489	66596738	chr18	73974071	73974320	chr18	56064244	56064493
chr18	66601231	66601480	chr18	73985549	73985798	chr18	56148666	56148915
chr18	69966208	69966457	chr18	73987174	73987423	chr18	56149037	56149286
chr18	70088939	70089188	chr18	48466086	48466335	chr18	56244171	56244420
chr18	70093364	70093613	chr18	48928528	48928777	chr18	56245242	56245491
chr18	70107935	70108184	chr18	48929078	48929327	chr18	56247063	56247312
chr18	70198868	70199117	chr18	48951524	48951773	chr18	56527049	56527298
chr18	70200128	70200377	chr18	48951856	48952105	chr18	56528011	56528260
chr18	70200907	70201156	chr18	48959882	48960131	chr18	56532190	56532439
chr18	70674976	70675225	chr18	48960304	48960553	chr18	56532882	56533131
chr18	70682967	70683216	chr18	49004516	49004765	chr18	56534208	56534457
chr18	73541191	73541440	chr18	49006425	49006674	chr18	56601487	56601736
chr18	73548315	73548564	chr18	49017695	49017944	chr18	56633131	56633380
chr18	73605288	73605537	chr18	49148010	49148259	chr18	56635259	56635508
chr18	73704485	73704734	chr18	49149627	49149876	chr18	56719296	56719545
chr18	73752546	73752795	chr18	49149933	49150182	chr18	58037828	58038077
chr18	73764755	73765004	chr18	49150459	49150708	chr18	58038328	58038577
chr18	73768420	73768669	chr18	50239403	50239652	chr18	58039111	58039360
chr18	73775027	73775276	chr18	50242789	50243038	chr18	58040290	58040539
chr18	73788751	73789000	chr18	50243069	50243318	chr18	59424085	59424334
chr18	73851108	73851357	chr18	50243367	50243616	chr18	59424765	59425014
chr18	73854713	73854962	chr18	50247229	50247478	chr18	59441324	59441573

FIGURE 2 (Page 11 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr18	59499718	59499967	chr18	60401084	60401333	chr18	60974839	60975088
chr18	59500226	59500475	chr18	60412199	60412448	chr18	60984378	60984627
chr18	59507920	59508169	chr18	60420053	60420302	chr18	61055997	61056246
chr18	59543044	59543293	chr18	60434561	60434810	chr18	61088669	61088918
chr18	59556453	59556702	chr18	60449047	60449296	chr18	61090462	61090711
chr18	59614976	59615225	chr18	60452471	60452720	chr18	61100325	61100574
chr18	59625553	59625802	chr18	60456757	60457006	chr18	61154459	61154708
chr18	59631090	59631339	chr18	60467147	60467396	chr18	61162901	61163150
chr18	59631779	59632028	chr18	60469757	60470006	chr18	61170629	61170878
chr18	59632388	59632637	chr18	60514168	60514417	chr18	61197492	61197741
chr18	59695718	59695967	chr18	60537888	60538137	chr18	61213452	61213701
chr18	59696087	59696336	chr18	60544380	60544629	chr18	61213900	61214149
chr18	59696563	59696812	chr18	60544742	60544991	chr18	61216749	61216998
chr18	59712419	59712668	chr18	60545407	60545656	chr18	61232928	61233177
chr18	59756979	59757228	chr18	60611934	60612183	chr18	61236079	61236328
chr18	59770789	59771038	chr18	60617081	60617330	chr18	61253277	61253526
chr18	59805993	59806242	chr18	60617641	60617890	chr18	61259478	61259727
chr18	59813876	59814125	chr18	60630438	60630687	chr18	61264533	61264782
chr18	59855411	59855660	chr18	60630849	60631098	chr18	61289440	61289689
chr18	59909054	59909303	chr18	60639840	60640089	chr18	61368889	61369138
chr18	59909336	59909585	chr18	60640300	60640549	chr18	61370663	61370912
chr18	59910140	59910389	chr18	60776944	60777193	chr18	63685629	63685878
chr18	59910602	59910851	chr18	60813738	60813987	chr18	63711583	63711832
chr18	59973803	59974052	chr18	60814585	60814834	chr18	63721183	63721432
chr18	59996728	59996977	chr18	60938666	60938915	chr18	64379625	64379874
chr18	60010541	60010790	chr18	60939202	60939451	chr18	64398368	64398617
chr18	60011381	60011630	chr18	60950541	60950790	chr18	64410092	64410341
chr18	60240694	60240943	chr18	60953225	60953474	chr18	64412606	64412855
chr18	60378353	60378602	chr18	60958474	60958723	chr18	64417267	64417516
chr18	60384593	60384842	chr18	60960879	60961128	chr18	64421562	64421811
chr18	60384979	60385228	chr18	60965583	60965832	chr18	64426871	64427120
chr18	60389544	60389793	chr18	60970051	60970300	chr18	64436287	64436536

FIGURE 2 (Page 12 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr18	64593727	64593976	chr18	25279010	25279259	chr18	25682079	25682328
chr18	24962594	24962843	chr18	25283817	25284066	chr18	25689906	25690155
chr18	24964748	24964997	chr18	25292042	25292291	chr18	25700540	25700789
chr18	25003553	25003802	chr18	25295429	25295678	chr18	25703648	25703897
chr18	25033345	25033594	chr18	25330345	25330594	chr18	25707188	25707437
chr18	25037434	25037683	chr18	25332963	25333212	chr18	25716100	25716349
chr18	25052992	25053241	chr18	25333284	25333533	chr18	25722951	25723200
chr18	25080042	25080291	chr18	25345252	25345501	chr18	25725055	25725304
chr18	25124915	25125164	chr18	25384702	25384951	chr18	25746105	25748354
chr18	25128057	25128306	chr18	25510767	25511016	chr18	26214441	26214690
chr18	25140257	25140506	chr18	25519108	25519357	chr18	26236788	26237037
chr18	25140722	25140971	chr18	25531724	25531973	chr18	26261336	26261585
chr18	25147759	25148008	chr18	25543168	25543417	chr18	26265937	26266186
chr18	25157828	25158077	chr18	25547425	25547674	chr18	26285286	26285535
chr18	25160244	25160493	chr18	25557343	25557592	chr18	26289006	26289255
chr18	25166354	25166603	chr18	25564483	25564732	chr18	26308405	26308654
chr18	25172986	25173235	chr18	25583087	25583336	chr18	26308962	26309211
chr18	25174972	25175221	chr18	25585083	25585332	chr18	26312054	26312303
chr18	25175640	25175889	chr18	25589902	25590151	chr18	26316149	26316398
chr18	25182855	25183104	chr18	25601966	25602215	chr18	47100905	47101154
chr18	25198531	25198780	chr18	25611986	25612235	chr18	47130278	47130527
chr18	25201886	25202135	chr18	25617498	25617747	chr18	47366936	47367185
chr18	25209395	25209644	chr18	25629105	25629354	chr18	47371131	47371380
chr18	25213880	25214129	chr18	25631180	25631429	chr18	47689868	47690117
chr18	25216383	25217137	chr18	25643639	25643888	chr18	47690267	47690516
chr18	25228617	25228866	chr18	25645371	25645620	chr18	47703689	47703938
chr18	25231535	25231784	chr18	25665025	25665274	chr18	47708139	47708388
chr18	25245662	25245911	chr18	25666614	25666863	chr18	47717445	47717694
chr18	25246300	25246549	chr18	25667057	25667306	chr18	47738531	47738780
chr18	25254843	25255092	chr18	25670741	25670990	chr18	47777438	47777687
chr18	25265476	25265725	chr18	25673208	25673457	chr18	47794212	47794461
chr18	25278532	25278781	chr18	25681018	25681267	chr18	47799400	47799649

FIGURE 2 (Page 13 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr18	47802662	47802911	chr18	43918519	43918768	chr18	46599847	46600096
chr18	39874986	39875235	chr18	43918998	43919247	chr18	46601866	46602115
chr18	39891490	39891739	chr18	43919544	43919793	chr18	46983278	46983527
chr18	42204240	42204489	chr18	44086123	44086372	chr18	47013103	47013352
chr18	42313146	42313395	chr18	44094299	44094548	chr18	47016802	47017051
chr18	42313939	42314188	chr18	44167214	44167463	chr18	286722	286971
chr18	42314302	42314551	chr18	44167791	44168040	chr18	298393	298642
chr18	42349923	42350172	chr18	44168296	44168545	chr18	338127	338376
chr18	42407165	42407414	chr18	44173512	44173761	chr18	340427	340676
chr18	42418033	42418282	chr18	44174626	44174875	chr18	387135	387384
chr18	42449702	42449951	chr18	44857765	44858014	chr18	391175	391424
chr18	42450400	42450649	chr18	45110525	45110774	chr18	434440	434689
chr18	42461427	42461676	chr18	451197621	451197870	chr18	457527	457776
chr18	42462834	42463083	chr18	45197871	45198120	chr18	460384	460633
chr18	42468595	42468844	chr18	45198760	45199009	chr18	469281	469530
chr18	42468938	42469187	chr18	45234886	45235135	chr18	572716	572965
chr18	42644003	42644252	chr18	45235523	45235772	chr18	598352	598601
chr18	42644494	42644743	chr18	45238370	45238619	chr18	615413	615662
chr18	42645387	42645636	chr18	45335896	45336145	chr18	618280	618529
chr18	42741059	42741308	chr18	45336353	45336602	chr18	901111	901360
chr18	42745733	42745982	chr18	45357461	45357710	chr18	939383	939632
chr18	42746460	42746709	chr18	45357745	45357994	chr18	26321283	26321532
chr18	42746911	42747160	chr18	45358502	45358751	chr18	26337865	26338114
chr18	42747473	42747722	chr18	45362379	45362628	chr18	26345717	26345966
chr18	43796659	43796908	chr18	45397588	45397837	chr18	26365387	26365636
chr18	43820471	43820720	chr18	45864183	45864432	chr18	26367518	26367767
chr18	43841913	43842162	chr18	46464382	46464631	chr18	26392500	26392749
chr18	43842202	43842451	chr18	46473428	46473677	chr18	26406499	26406748
chr18	43842647	43842896	chr18	46474545	46474794	chr18	26711590	26711839
chr18	43845701	43845950	chr18	46481797	46482046	chr18	26729342	26729591
chr18	43846353	43846602	chr18	46483016	46483265	chr18	26733981	26734230
chr18	43917652	43917901	chr18	46589014	46589263	chr18	26739336	26739585

FIGURE 2 (Page 14 of 31)

chromosome	start	stop
chr18	26750685	26750934
chr18	26752487	26752736
chr18	26757573	26757822
chr18	27244230	27244479
chr18	35076347	35076596
chr18	35147912	35148161
chr18	35282655	35282904
chr18	35292141	35292390
chr18	35304247	35304496
chr18	35311302	35311551
chr18	36463709	36463958
chr18	36473033	36473282
chr18	36475403	36475652
chr18	36479145	36479394
chr18	36482038	36482287
chr18	36483959	36484208
chr18	36485761	36486010
chr18	36487627	36487876
chr18	36489135	36489384
chr18	36531244	36531493
chr18	36554006	36554255
chr18	36558050	36558299
chr18	36569619	36569868
chr18	1017928	1018177
chr18	241538	241787
chr18	242831	243080

FIGURE 2 (Page 15 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr13	19520254	19520499	chr13	26939424	26939672	chr13	32759010	32759261
chr13	19523790	19524042	chr13	26939942	26940196	chr13	32784101	32784351
chr13	20531155	20531403	chr13	26952790	26953038	chr13	32784570	32784820
chr13	20542998	20543246	chr13	26956549	26956797	chr13	32786822	32787069
chr13	21182580	21182828	chr13	26957204	26957453	chr13	32787129	32787381
chr13	22255353	22255601	chr13	26957736	26957984	chr13	32811213	32811463
chr13	26795002	26795250	chr13	26971501	26971749	chr13	32931779	32932029
chr13	26795380	26795628	chr13	26974853	26975100	chr13	33016510	33016761
chr13	26809864	26810112	chr13	26976810	26977059	chr13	33016929	33017182
chr13	26810422	26810669	chr13	27117460	27117707	chr13	33017432	33017681
chr13	26810951	26811198	chr13	27216274	27216526	chr13	36342494	36342743
chr13	26811424	26811674	chr13	29111011	29111260	chr13	36343935	36344179
chr13	26820545	26820794	chr13	29233772	29234021	chr13	36344330	36344577
chr13	26824800	26825048	chr13	29310139	29310388	chr13	36345572	36345822
chr13	26825155	26825405	chr13	29563140	29563389	chr13	36346001	36346253
chr13	26827150	26827400	chr13	31499212	31499458	chr13	36349204	36349453
chr13	26827690	26827939	chr13	31499519	31499772	chr13	36354263	36354513
chr13	26829149	26829396	chr13	31609109	31609361	chr13	36358322	36358572
chr13	26829663	26829911	chr13	31653571	31653822	chr13	36361603	36361850
chr13	26830241	26830491	chr13	32149129	32149380	chr13	36362050	36362297
chr13	26830739	26830987	chr13	32169142	32169392	chr13	36362531	36362778
chr13	26834633	26834880	chr13	32173373	32173622	chr13	36363439	36363689
chr13	26836514	26836761	chr13	32271292	32271539	chr13	36366103	36366353
chr13	26843066	26843315	chr13	32316109	32316362	chr13	36387184	36387431
chr13	26851259	26851507	chr13	32316651	32316898	chr13	39584383	39584630
chr13	26864367	26864615	chr13	32339360	32339608	chr13	39585052	39585299
chr13	26888243	26888485	chr13	32359633	32359883	chr13	39585675	39585922
chr13	26888660	26888907	chr13	32433382	32433629	chr13	40090054	40090298
chr13	26898691	26898936	chr13	32434990	32435237	chr13	40093695	40093942
chr13	26925436	26925684	chr13	32600409	32600659	chr13	40098460	40098707
chr13	26926976	26927224	chr13	32653613	32653857	chr13	40098913	40099157
chr13	26938928	26939175	chr13	32758635	32758877	chr13	40101053	40101303

FIGURE 2 (Page 16 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr13	40101582	40101832	chr13	68311152	68311402	chr13	73320095	73320399
chr13	40384942	40385189	chr13	68316131	68316378	chr13	77198180	77198433
chr13	40391440	40391687	chr13	68428699	68428952	chr13	77422635	77422880
chr13	40392221	40392468	chr13	68429013	68429258	chr13	77449372	77449622
chr13	40398309	40398557	chr13	68570544	68570791	chr13	77454212	77454458
chr13	40398589	40398840	chr13	71269531	71269778	chr13	77454985	77455228
chr13	40406460	40406711	chr13	71454173	71454420	chr13	77455721	77455971
chr13	40420794	40421037	chr13	71456180	71456429	chr13	77456394	77456643
chr13	40613790	40614037	chr13	71463069	71463321	chr13	77476090	77476339
chr13	42846940	42847187	chr13	71464674	71464887	chr13	77477580	77477833
chr13	42876485	42876732	chr13	71465060	71465307	chr13	77513439	77513692
chr13	42896653	42896897	chr13	71466315	71466560	chr13	77591495	77591742
chr13	43149529	43149781	chr13	71466724	71466967	chr13	77595913	77596160
chr13	43164772	43165022	chr13	71479089	71479339	chr13	77598010	77598257
chr13	43177841	43178091	chr13	71479359	71479608	chr13	77630724	77630968
chr13	52163649	52163899	chr13	71619104	71619349	chr13	77633431	77633680
chr13	52541402	52541652	chr13	71750165	71750409	chr13	77636364	77636611
chr13	52549279	52549529	chr13	71750461	71750708	chr13	77639621	77639873
chr13	53323273	53323520	chr13	71762243	71762493	chr13	77642789	77643041
chr13	53350485	53350727	chr13	71762640	71762893	chr13	77652265	77652509
chr13	53505942	53506191	chr13	71781780	71782031	chr13	77652705	77652952
chr13	53555869	53556121	chr13	71783224	71783469	chr13	77675169	77675423
chr13	53606740	53606987	chr13	71784572	71784822	chr13	77699605	77699849
chr13	53610142	53610389	chr13	73226435	73226679	chr13	77713099	77713352
chr13	53619813	53620063	chr13	73233924	73234168	chr13	77714235	77714481
chr13	53622959	53623209	chr13	73251601	73251851	chr13	77742675	77742920
chr13	59239399	59239649	chr13	73261064	73261309	chr13	77759281	77759528
chr13	62064603	62064852	chr13	73284370	73284617	chr13	78540689	78540943
chr13	68236363	68236610	chr13	73296842	73297089	chr13	89222075	89222318
chr13	68238595	68238842	chr13	73305479	73305731	chr13	94063232	94063481
chr13	68241140	68241393	chr13	73305843	73306093			
chr13	68310701	68310948	chr13	73309111	73309357			

FIGURE 2 (Page 17 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr18	11701840	11702094	chr18	13040411	13040657	chr18	29705097	29705347
chr18	11742882	11743124	chr18	13049051	13049304	chr18	29707000	29707248
chr18	11743713	11743964	chr18	13053013	13053260	chr18	29711514	29711758
chr18	11743974	11744216	chr18	13057128	13057370	chr18	29772771	29773020
chr18	11745797	11746044	chr18	13725324	13725566	chr18	29773038	29773291
chr18	11766505	11766752	chr18	13731669	13731919	chr18	29790110	29790357
chr18	11784594	11784841	chr18	13930816	13931068	chr18	29794196	29794446
chr18	11784857	11785107	chr18	13936510	13936761	chr18	29814455	29814703
chr18	11789414	11789659	chr18	13939030	13939277	chr18	29816677	29816924
chr18	11791405	11791649	chr18	13941500	13941742	chr18	29824628	29824872
chr18	11831898	11832147	chr18	18709626	18709868	chr18	29941738	29941986
chr18	11840071	11840318	chr18	19345013	19345263	chr18	29942404	29942657
chr18	11853018	11853271	chr18	22314388	22314639	chr18	29947043	29947286
chr18	12786330	12786572	chr18	22317204	22317445	chr18	29963606	29963859
chr18	12793790	12794034	chr18	22318584	22318826	chr18	29967321	29967565
chr18	12813342	12813590	chr18	22474179	22474425	chr18	29988241	29988491
chr18	12813796	12814046	chr18	22503094	22503336	chr18	29989060	29989304
chr18	12836508	12836756	chr18	22514422	22514675	chr18	33702205	33702458
chr18	12837826	12838071	chr18	22527670	22527920	chr18	33703293	33703541
chr18	12839370	12839619	chr18	22713619	22713863	chr18	33703805	33704052
chr18	12839810	12840054	chr18	22715092	22715939	chr18	33710269	33710511
chr18	12840789	12841033	chr18	22746015	22746258	chr18	33831399	33831647
chr18	12955546	12955788	chr18	22751334	22751581	chr18	3384883	33849083
chr18	12968393	12968643	chr18	26362055	26362304	chr18	33955911	33956161
chr18	12984202	12984449	chr18	27204525	27205070	chr18	33967910	33968160
chr18	12984953	12985198	chr18	27683185	27683429	chr18	35703018	35703269
chr18	12986018	12986263	chr18	28529366	28529610	chr18	35705739	35705986
chr18	13007427	13007681	chr18	28538550	28538803	chr18	36487627	36487876
chr18	13007815	13008064	chr18	28550611	28550855	chr18	52428207	52428460
chr18	13019314	13019562	chr18	28551683	28551928	chr18	52472800	52473050
chr18	13029804	13030046	chr18	28584849	28585102	chr18	52702246	52702490
chr18	13036834	13037081	chr18	29704799	29705048	chr18	52703207	52703449

FIGURE 2 (Page 18 of 31)

chromosome	start	stop
chr18	52719979	52720226
chr18	53843255	53843507
chr18	53843510	53843764
chr18	53843879	53844121
chr18	53870090	53870332
chr18	53873286	53873528
chr18	53913040	53913287
chr18	54141132	54141379
chr18	54154870	54155117
chr18	54183137	54183381
chr18	54271758	54272011
chr18	54352988	54353238
chr18	60378353	60378602
chr18	61164144	61164393

chromosome	start	stop
1	0	223,700,000
2	0	223,700,000
3	0	223,700,000
4	0	223,700,000
5	0	223,700,000
6	0	223,700,000
7	0	223,700,000
8	0	223,700,000
9	0	223,700,000
10	0	223,700,000
11	0	223,700,000
12	0	223,700,000
13	0	223,700,000
14	0	223,700,000
15	0	223,700,000
16	0	223,700,000
17	0	223,700,000
18	0	223,700,000
19	0	223,700,000
20	0	223,700,000
21	0	223,700,000
22	0	223,700,000
X	0	223,700,000
Y	0	223,700,000

chromosome	start	stop
------------	-------	------

FIGURE 2 (Page 19 of 31)

chromosome	start	stop	chromosome	start	stop	chromosome	start	stop
chr21	15774354	15774598	chr21	17160410	17160656	chr21	30568450	30568700
chr21	15774701	15774951	chr21	17160972	17161216	chr21	30901196	30901442
chr21	15941894	15942144	chr21	19300496	19300748	chr21	30950156	30950403
chr21	15942369	15942621	chr21	19396568	19396818	chr21	30950626	30950878
chr21	16330536	16330783	chr21	25701930	25702183	chr21	30952977	30953227
chr21	16331481	16331728	chr21	25702313	25702560	chr21	30953281	30953523
chr21	16332238	16332485	chr21	27985131	27985378	chr21	30953836	30954083
chr21	16333694	16333941	chr21	27989556	27989803	chr21	30966825	30967077
chr21	16335023	16335271	chr21	28000089	28000336	chr21	30968330	30968574
chr21	16336651	16336901	chr21	28007209	28007456	chr21	30969659	30969903
chr21	16337562	16337816	chr21	28007550	28007793	chr21	31015490	31015733
chr21	16338381	16338633	chr21	28008308	28008555	chr21	31022236	31022488
chr21	16338902	16339154	chr21	28054591	28054835	chr21	31091818	31092065
chr21	16339772	16340019	chr21	28292245	28292493	chr21	31109787	31110039
chr21	16340409	16340656	chr21	28294400	28294644	chr21	31110149	31110396
chr21	16344480	16344730	chr21	28308980	28309224	chr21	31113996	31114249
chr21	16349101	16349354	chr21	28311566	28311813	chr21	31120815	31121064
chr21	16362700	16362954	chr21	28316801	28317045	chr21	31144588	31144835
chr21	16375294	16375539	chr21	29796060	29796309	chr21	35154438	35154687
chr21	16376059	16376308	chr21	30112410	30112660	chr21	40549354	40549598
chr21	16376695	16376942	chr21	30112772	30113024	chr21	40759240	40759493
chr21	16519760	16520004	chr21	30113106	30113350	chr21	40839325	40839569
chr21	16529308	16529552	chr21	30114870	30115119	chr21	40846844	40847093
chr21	16529940	16530189	chr21	30301073	30301322	chr21	41160090	41160343
chr21	16531755	16532002	chr21	30301466	30301713	chr21	41171755	41171999
chr21	16533111	16533355	chr21	30317970	30318220	chr21	41298214	41298467
chr21	16537385	16537631	chr21	30358601	30358851	chr21	41592339	41592583
chr21	16545527	16545774	chr21	30397555	30397801	chr21	41596188	41596432
chr21	16798480	16798722	chr21	30419470	30419720	chr21	41741900	41742144
chr21	16792457	16792708	chr21	30527461	30527714	chr21	43219112	43219356
chr21	17149630	17149877	chr21	30528513	30528763	chr21	46468654	46468898
chr21	45114956	45115203	chr21	30554566	30554808	chr21	46495701	46495954

FIGURE 2 (Page 20 of 31)




FIGURE 2 (Page 21 of 31)

chromosome	start	stop
chr2	50176039	50176288
chr2	63596895	63597144
chr2	103538910	103539159

chromosome	start	stop
------------	-------	------

chromosome	start	stop
------------	-------	------

FIGURE 2 (Page 22 of 31)




FIGURE 2 (Page 23 of 31)

chromosome	start	stop
chr4	21947614	21947863
chr4	62556697	62556946
chr4	72418734	72418983
chr4	122570216	122570455
chr4	130264593	130264842

chromosome	start	stop
------------	-------	------

chromosome	start	stop
------------	-------	------

FIGURE 2 (Page 24 of 31)

chromosome	start	stop
chr5	59283789	59284038
chr5	88597843	88598092
chr5	163891445	163891694
chr5	163896919	163897168

chromosome	start	stop
------------	-------	------

chromosome	start	stop
1	1000000	1000000
2	1000000	1000000
3	1000000	1000000
4	1000000	1000000
5	1000000	1000000
6	1000000	1000000
7	1000000	1000000
8	1000000	1000000
9	1000000	1000000
10	1000000	1000000
11	1000000	1000000
12	1000000	1000000
13	1000000	1000000
14	1000000	1000000
15	1000000	1000000
16	1000000	1000000
17	1000000	1000000
18	1000000	1000000
19	1000000	1000000
20	1000000	1000000
21	1000000	1000000
22	1000000	1000000
23	1000000	1000000
24	1000000	1000000
25	1000000	1000000
26	1000000	1000000
27	1000000	1000000
28	1000000	1000000
29	1000000	1000000
30	1000000	1000000
31	1000000	1000000
32	1000000	1000000
33	1000000	1000000
34	1000000	1000000
35	1000000	1000000
36	1000000	1000000
37	1000000	1000000
38	1000000	1000000
39	1000000	1000000
40	1000000	1000000
41	1000000	1000000
42	1000000	1000000
43	1000000	1000000
44	1000000	1000000
45	1000000	1000000
46	1000000	1000000
47	1000000	1000000
48	1000000	1000000
49	1000000	1000000
50	1000000	1000000
51	1000000	1000000
52	1000000	1000000
53	1000000	1000000
54	1000000	1000000
55	1000000	1000000
56	1000000	1000000
57	1000000	1000000
58	1000000	1000000
59	1000000	1000000
60	1000000	1000000
61	1000000	1000000
62	1000000	1000000
63	1000000	1000000
64	1000000	1000000
65	1000000	1000000
66	1000000	1000000
67	1000000	1000000
68	1000000	1000000
69	1000000	1000000
70	1000000	1000000
71	1000000	1000000
72	1000000	1000000
73	1000000	1000000
74	1000000	1000000
75	1000000	1000000
76	1000000	1000000
77	1000000	1000000
78	1000000	1000000
79	1000000	1000000
80	1000000	1000000
81	1000000	1000000
82	1000000	1000000
83	1000000	1000000
84	1000000	1000000
85	1000000	1000000
86	1000000	1000000
87	1000000	1000000
88	1000000	1000000
89	1000000	1000000
90	1000000	1000000
91	1000000	1000000
92	1000000	1000000
93	1000000	1000000
94	1000000	1000000
95	1000000	1000000
96	1000000	1000000
97	1000000	1000000
98	1000000	1000000
99	1000000	1000000
100	1000000	1000000

FIGURE 2 (Page 25 of 31)

chromosome	start	stop
chr6	8529895	8530144
chr6	96221217	96221466
chr6	112441026	112441275
chr6	118481182	118481431
chr6	129800139	129800388

FIGURE 2 (Page 26 of 31)

chromosome	start	stop
chr7	22725339	22725568
chr7	42250786	42251035

FIGURE 2 (Page 27 of 31)

chromosome	start	stop
chr9	2600491	2600740
chr9	14299731	14299980
chr9	72782230	72782479
chr9	82880926	82881175

FIGURE 2 (Page 28 of 31)

chromosome	start	stop
chr10	52560407	52560656

chromosome	start	stop

chromosome	start	stop

FIGURE 2 (Page 29 of 31)

chromosome	start	stop
chr11	10786015	10786264

chromosome	start	stop

chromosome	start	stop


FIGURE 2 (Page 30 of 31)

chromosome	start	stop
chr12	29714971	29715220
chr12	40732964	40733213

chromosome	start	stop

chromosome	start	stop

FIGURE 2 (Page 31 of 31)



Sample
FIGURE 3







FIGURE 4




Sample
FIGURE 5




Sample
FIGURE 6



Sample
FIGURE 7





Sample
FIGURE 9




FIGURE 10

