



(51) International Patent Classification:
G06F 19/00 (2011.01)

(21) International Application Number:

PCT/US20 11/042416

(22) International Filing Date:

29 June 2011 (29.06.2011)

(25) Filing Language:

English

(26) **Publication Language:**

English

(30) Priority Data:

61/359.708

29 June 2010 (29.06.2010)

II S

(71) Applicant (for all designated States except US): **THE UNIVERSITY OF NORTH CAROLINA AT CHAPEL HILL** [US/US]; 308 Bynum Hall, Campus Box 4105, Chapel Hill, N.C. 27599-4105 (US).

(72) Inventors; and

(75) **Investors/ Applicants** (for US only): **KOCIS, Keith, Christopher** [US/US]; 507 Highgrove Drive, Chapel Hill, N C 275 16 (US). **KOCIS, Daniel, Joseph** [US/US]; 12 Harmony Street, Manorville, NY 11949 (US).

(74) **Agent: HUNT, Gregory, A.;** Jenkins, Wilson, Taylor & Hunt, P.A., Suite 1200, University Tower, 3100 Tower Boulevard, Durham, N C 27707 (US).


(81) Designated States (*unless otherwise indicated, for every kind of national protection available*): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IS, JP, KE, KG, KM, KN, KP, KR, KZ, LA, LC, LK, LR, LS, LT, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PE, PG, PH, PL, PT, RO, RS, RU, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.

(84) Designated States (*unless otherwise indicated, for every kind of regional protection available*): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

Published:

- without international search report and to be republished upon receipt of that report (Rule 48.2(q))

(54) Title: METHODS, SYSTEMS, AND COMPUTER READABLE MEDIA FOR EVALUATING A HOSPITAL PATIENT'S RISK OF MORTALITY



ETG I

(57) Abstract: The subject matter includes methods, systems, and computer readable media for evaluating a hospital patient's risk of mortality. According to one aspect of the subject matter described herein, a method for evaluating a hospital patient's risk of mortality is provided. The method includes collecting data from physiologic signals generated by patient monitors, physiologic signals of organ function, and demographic information for a patient. The method further includes determining a measure of the variability of at least one of the physiologic signals. The method further includes analyzing data and the measure of variability to determine whether a value for a particular physiologic or demographic variable falls within a critical interval for the variable that indicates that the value is predictive of mortality or survival. Each time a value for a physiological or demographic variable for the patient falls within a critical interval, the method includes recording the occurrence of an event for the patient. The method further includes counting the number of events for the patient over a time period. The method further includes, generating, based on the account output perceptible by human user that indicates the patient's risk of mortality or likelihood of survival.

DESCRIPTION

METHODS, SYSTEMS, AND COMPUTER READABLE MEDIA FOR
EVALUATING A HOSPITAL PATIENTS RISK OF MORTALITY

PRIORITY CLAIM

This application claims the benefit of U.S. Provisional Patent Application
5 Serial No. 61/359,708, filed June 29, 2010; the disclosure of which is
incorporated herein by reference in its entirety.

GOVERNMENT INTEREST

This presently disclosed subject matter was made with U.S. Government
10 support under Grant No. 1R41HD049935-01 awarded by National Institutes of
Health. Thus, the U.S. Government has certain rights in the presently
disclosed subject matter.

TECHNICAL FIELD

15 The subject matter described herein relates to analyzing data regarding
hospital patients. More particularly, the subject matter described herein relates
to methods, systems, and computer readable media for evaluating a hospital
patient's risk of mortality.

20 BACKGROUND

By the year 2010 there will be approximately 5,500 pediatric intensive
care unit (PICU) beds, about 22,000 neonatal intensive care unit (NICU) beds,
and about 110,000 non-intensive care (non-ICU) pediatric care beds in the
United States of America. Rendering care to these hospitalized children
25 represents a large, significant, and growing cost to health care expense. The
importance of intensive care and non-intensive care monitoring in patient safety
is recognized as a priority by the Joint Commission on Accreditation of Hospital
Organizations. Intensive care and non-intensive care monitors are used to
continuously evaluate the clinical status of patients and track their response to
30 a wide range of interventions. Despite tremendous advancements in computer
technology and bioinformatics, intensive care and non-intensive care
monitoring devices have not provided significantly new information to bedside

caregivers. Simplified cardiopulmonary monitoring is now routinely performed in hospitalized non-intensive care patients.

Rapid Response (RR) teams have been established to rescue non-ICU pediatric patients who are decompensating and in need of critical care evaluation. These teams are triggered by any health care provider (nurse, respiratory therapist, physician, etc.) or family members resulting in a stat evaluation by an expert team of health care providers that can initiate lifesaving support and transfer patients to an ICU. These teams have been shown to be effective at saving the lives of both pediatric and adult hospitalized inpatients.

Software tools are available for evaluating a patient's risk of mortality for patients in PICU, NICU and non-ICU pediatric facilities. However, conventional tools only evaluate a patient's risk of mortality based on raw variables considered individually. Models that evaluate a patient's risk of mortality based on raw variables considered individually lack accuracy and lead to false alarms.

For example, conventional monitors that monitor intensive care patients' vital signs alarm frequently due to variability in the raw data, even when patients are not at risk of imminent death. As a result, hospital staff and physicians often ignore such alarms.

Accordingly, there exists a long felt need for methods, systems, and computer readable media for evaluating a hospital patient's risk of mortality.

SUMMARY

The subject matter includes methods, systems, and computer readable media for evaluating a hospital patient's risk of mortality. According to one aspect of the subject matter described herein, a method for evaluating a hospital patient's risk of mortality is provided. The method includes collecting data from physiologic signals generated by patient monitors, physiologic signals of organ function, and demographic information for a patient. The method further includes determining a measure of the variability of at least one of the physiologic signals. The method further includes analyzing data and the measure of variability to determine whether a value for a particular physiologic or demographic variable falls within a critical interval for the variable that indicates that the value is predictive of mortality or survival. Each time a value

for a physiological or demographic variable for the patient falls within a critical interval, the method includes recording the occurrence of an event for the patient. The method further includes counting the number of events for the patient over a time period. The method further includes, generating, based on
5 the account output perceptible by human user that indicates the patient's risk of mortality or likelihood of survival.

The subject matter described herein can be implemented in software in combination with hardware and/or firmware. For example, the subject matter described herein can be implemented in software executed by a processor. In
10 one exemplary implementation, the subject matter described herein can be implemented using a non-transitory computer readable medium having stored thereon computer executable instructions that when executed by the processor of a computer control the computer to perform steps. Exemplary computer readable media suitable for implementing the subject matter described herein
15 include non-transitory computer-readable media, such as disk memory devices, chip memory devices, programmable logic devices, and application specific integrated circuits. In addition, a computer readable medium that implements the subject matter described herein may be located on a single device or computing platform or may be distributed across multiple devices or computing
20 platforms.

BRIEF DESCRIPTION OF THE DRAWINGS

Preferred embodiments of the subject matter described herein will now be explained with reference to the accompanying drawings, wherein like
25 reference numerals represent like parts, of which:

Figure 1 is a block diagram illustrating exemplary overall components of an embodiment of the subject matter described herein;

Figure 2 is a block diagram illustrating exemplary components of a system for evaluating a hospital patient's risk of mortality according to an
30 embodiment of the subject matter described herein;

Figure 3 is a flow chart illustrating exemplary steps for evaluating a hospital patient's risk of mortality according to an embodiment of the subject matter described herein;

Figures 4A and 4B are graphs of ECG voltage versus time that may be analyzed by the system illustrated in Figure 2 to evaluate a patient's risk of mortality according to an embodiment of the subject matter described herein;

Figures 5A and 5B are tables respectively illustrating statistical measures for variables for live and dead PICU patients that may be used by the system illustrated in Figure 2 to evaluate a given patient's risk of mortality according to an embodiment of the subject matter described herein;

Figures 6A and 6B are tables illustrating statistical measures for additional variables for dead PICU and live PICU patients that may be used by the system illustrated in Figure 2 to evaluate a given patient's risk of mortality according to an embodiment of the subject matter described herein;

Figure 7 is a table illustrating the dividing of data points into critical intervals for a PICU patient according to an embodiment of the subject matter described herein;

Figures 8A and 8B are tables illustrating exemplary statistical measures for different variables for dead PICU and live PICU patients according to an embodiment of the subject matter described herein;

Figures 9A and 9B are tables illustrating exemplary measures of patient physiological data that may be obtained and used by the system illustrated in Figure 2 to evaluate a given patient's risk of mortality according to an embodiment of the subject matter described herein;

Figures 10A and 10B illustrate variables, intervals, and statistical measures that may be used by the system illustrated in Figure 2 to evaluate a hospital patient's risk of mortality according to an embodiment of the subject matter described herein;

Figures 11A and 11B illustrate intervals that may be used by the system in Figure 2 to evaluate a patient's risk of mortality according to an embodiment of the subject matter described herein;

Figure 12A is a table, and Figure 12B is a graph illustrating critical intervals for live and dead patients for a particular variable that may be used by the system illustrated in Figure 2 to evaluate a patient's risk of mortality according to an embodiment of the subject matter described herein;

Figures 13A and 13B are plots of sodium frequency versus time for a given patient according to an embodiment of the subject matter described herein;

Figures 14A-14C are tables illustrating counts of multi-feed
5 transformations on expanded PICU data where derivatives and second derivatives are computed to determine measures of variability of the variable being analyzed;

Figure 15A is a table, and Figure 15B is a graph illustrating critical intervals of a particular variable that may be used by the system illustrated in
10 Figure 2 to evaluate a patient's risk of mortality according to an embodiment of the subject matter described herein;

Figures 15C is a table and Figure 15D is a graph of critical intervals for a particular variable that may be used by the system illustrated in Figure 2 to evaluate a patient's risk of mortality according to an embodiment of the subject
15 matter described herein;

Figure 16A is a graph of intermittent raw data;

Figure 16B is a graph of extrapolated data in 2 minute increments which calibrates all data into exact 2minute intervals; and

Figure 17A is a graph of patient event counts per organ and Figure 17B
20 is a table of the same data that may be generated by the system illustrated in Figure 2 according to an embodiment of the subject matter described herein;

Figure 17C includes graphs tracking each two minute score along with organ specific critical interval counts that may be generated by the system illustrated in Figure 2 according to an embodiment of the subject matter
25 described herein;

Figure 18A is a partial predictive model markup language (PMML) listing of a model for evaluating risk of mortality for a PICU patient according to an embodiment of the subject matter described herein;

Figure 18B is a partial PMML listing of a model for evaluating risk of
30 mortality for a rapid response (RR) patient according to an embodiment of the subject matter described herein;

Figure 19 is an example of an ECG signal with corruption that can be corrected using covariance according to an embodiment of the subject matter described herein;

Figure 20 is a graph of a QRS template;

5 Figure 21 is a graph of the covariance between the QRS template illustrated in Figure 20 and the contaminated ECG signal illustrated in Figure 19;

Figure 22 is a graph of the output of a quadrature-mirrored, second order linear infinite impulse response filter use to remove baseline shifts in a QRS wave according to an embodiment of the subject matter described herein;

Figure 23 is a graph of the R wave determined using a 2 sigma statistical filter according to an embodiment of the subject matter described herein; and

Figure 24 is a graph of the output of a threshold and first derivative filter, which represents the ECG signal with contamination removed and after post-processing according to an embodiment of the subject matter described herein.

DETAILED DESCRIPTION

In accordance with the subject matter disclosed herein, systems, methods, and computer readable media are provided for evaluating a hospital patient's risk of mortality. One exemplary implementation of the subject matter described herein is referred to as REALTROMINS, which is an acronym for REAL Time Risk of Mortality and INStability. REALTROMINS is a medical device for critically ill patients that integrates and analyzes numerous inputs: 1) real time continuous physiologic signals (e.g. electrocardiogram); 2) advanced measures of variability of these signals (e.g. spectral analysis); 3) physiologic based measures of organ function (e.g. serum glucose) and their variability; and 4) demographic and diagnosis related predictors of mortality; creating a real time, continuously updated risk of mortality score to aid clinicians in administering medical care. Described herein are REALTROMINS devices and methods for use in patient populations comprising premature infants in a neonatal intensive care unit (NICU; Neonatal REALTROMINS) and hospitalized children outside the ICU but in a non-ICU pediatric care facility and/or under the care and supervision of a Rapid Response (RR) team (Rapid Response

REALTROMINS). Also described herein are REALTROMINS devices and methods for use in the PICU. The REALTROMINS devices and methods can guide medical decision making resulting in earlier initiation or withdrawal of costly and risky therapeutic interventions to improve patient outcomes. In addition, the presently disclosed subject matter can decrease healthcare costs by better matching expensive and limited: 1) human resources (physician, nursing, etc.) and 2) hospital facilities (intensive care beds) to the changing needs of critically ill and hospitalized patients.

REALTROMINS is a bioinformatics technology designed to continuously assess the changing risk of mortality of critically ill pediatric patients and track the success or failure of ongoing medical interventions. Described herein are REALTROMINS technologies to be utilized in treating patients in PICU, NICU and non-ICU pediatric settings. Bedside caregivers can determine which therapies are leading to a decrease in risk while an increased mortality risk would trigger new approaches to the clinical problem. Patient outcomes can improve by displaying the overall risk score and the individual components contributing to that risk. By way of example and not limitation, heart rate variability has been shown to be highly predictive of mortality in a wide variety of critical illnesses in both adults and children. While not intended to be bound by any single theory or mechanism of action, the concept is that healthy biological systems are in a natural state of chaos and when critical illness develops, this natural variability is lost. Analysis of the low frequency (LF) and high frequency (HF) spectral content of heart rate variability, with calculation of a ratio of these two areas, has been predictive of mortality. As such, in some embodiments REALTROMINS couples these metrics, along with others, into an advanced multivariate analysis using SAS Enterprise Miner Software (SAS Institute, Cary, North Carolina, United States of America). In other embodiments, these metrics are assigned to critical intervals and used to generate hospital patient scores without using SAS Enterprise Miner Software.

Research Design and Methods

Projects and experiments were designed for the generation and collection of data for the development of the Neonatal REALTROMINS and

Rapid Response REALTROMINS. A prospective cohort study was used for the estimation and validation data sets. A sixteen-bed PICU (REALTROMINS), 45-bed NICU (Neonatal REALTROMINS), and a 69 non-ICU bed (Rapid Response REALTROMINS) tertiary care university based children's hospital was used in the studies. Annual admissions to the PICU, NICU, and non-ICU beds were approximately 1,000, 750, and 3,100, respectively. Patients ranged in age from 0 (as young as 22 weeks gestation) to 18 years, encompassing the full range of pediatric medical and surgical diagnoses. In the PICU, the male to female ratio is 1.27:1 and a diverse racial background is represented. The racial and ethnic breakdowns are as follows: 48% white, 31% black and 21% Hispanic and other races. Patient Enrollment Criteria: For the estimation (50% of patients) and validation (50% of patients) data sets, all patients admitted to the PICU, NICU, and non-ICU beds were enrolled over 9 months and comparison between survivors and non-survivors was examined. Since few patients died in the non-ICU beds additional outcome measures, such as transfer to the PICU, triggering a code blue (cardiopulmonary arrest) or Rapid Response team dispatch, were be examined.

Summary of Data Sources and Information Flow

As depicted in Figure 1, the following continuous physiologic signals were obtained in all subjects as determined by the clinical needs of the patient: 1) cardiac a) electrocardiogram; b) invasive arterial blood pressure; c) invasive central venous pressure; 2) respiratory a) respiratory plethysmograph; b) pulse oximetry; c) end tidal carbon dioxide concentrations; and 3) continuous temperature. These continuous digital data streams were written to a central server in an SQL database structure every 8 seconds. This database contained information detailing which monitor (i.e. patient) data is streaming from, UTC time stamps, and the individual data elements. All data were stored and backed up using standard procedures (RAID 5 hard drive configuration and digital tape).

The electrocardiogram (ECG) was digitized and a java based algorithm was used to detect the fiducial point (R wave of the ECG) from which the RR interval (time between R peaks), heart period, and heart rate (HR) time series

were derived. The front end loader manages data capture. This component consists of the open source Mirth HL-7 engine for capture of all HL-7 data feeds and a custom Java application that captures real-time ICU monitor data every 30 seconds (configurable) from the vendor's transactional database via a JDBC database connection. At the point of capture, the UTC time stamp, patient location, and patient IDs are standardized and data is then forwarded to the messaging system. A Java-based real-time patient system (RTPS) captures all data packets from the messaging system for pre-selected patient locations. The RTPS dynamically links all data packets for each patient and stores the data in-memory in a custom multivariate time series representation. The custom representation allows each variable to be stored at the variable's temporal resolution in addition to providing access to packets that provide all variable values during a fixed time window (typically 2 minutes). Signal processing and preprocessing algorithms are applied to the time series, which results in the creation of new time series variables. Decision support algorithms can then be applied on a continuous basis to the incoming time series data. Since a patient can be monitored for days to months at a time, the in-memory data is limited to the maximum time needed for signal processing and decision support while older data is pushed to a transactional database within the RTPS. A second Java application monitors the RTPS to detect patients that have been discharged from a location. Once detected, all data on that patient is extracted from the RTPS and moved to the data archival system. As part of this process, the patient data is converted from the custom time series representation to a format amendable to model development.

Spectral analyses were then performed. The HR time series was interpolated to achieve a uniform sampling interval and normalized to mean zero and variation from the mean. The amplitude spectrum was generated using the Fast Fourier transform (FFT). The power spectrum was derived from the amplitude spectra and used to determine frequency band specific (ultra-low, very low, low, and high frequency) variance of the HR time series for statistical tests of significance. Sequential records of Power spectra were generated, one every 128 beats (approximately every 1-2 minutes in children).

These results were then entered into the statistical and neural network software.

Physiologic based indices of organ function were obtained at an hourly frequency (i.e. urine output) or intermittently (i.e. serum glucose) according to patient care needs. Physiologic based measures of organ function included: 1) cardiac indicia, comprising non-invasive blood pressure (hourly systolic and diastolic), lactate level; 2) respiratory indicia, comprising respiratory rate, pH, **PCO₂**, PO₂, HCO₃, O₂ saturation, fractional inspired oxygen concentration (hourly); 3) neurologic indicia, comprising Glasgow coma score (hourly), pupillary reactions (hourly); 4) fluid indicia, comprising sodium, potassium, glucose, urea nitrogen, creatinine concentrations, ionized calcium, whole blood glucose measurements, urine output (hourly), hourly fluids administered (IV and/or enteral), weight on admission; 5) hematologic indicia, comprising hemoglobin, platelet count, prothrombin time-INR, partial thromboplastin time, D-dimers; 6) hepatic indicia, comprising bilirubin, transaminase, albumin, total protein; and 7) immunologic indicia, comprising temperature (hourly if not continuous), white cell count (total), absolute neutrophil count, absolute lymphocyte count, c reactive protein. Changes in these variables over time were investigated using the SAS procedure EXPAND to calculate first and second derivatives.

Demographic and diagnosis related predictors. 1) age; 2) general diagnostic groups (e.g., oncologic disease); 3) pre ICU risk factors (e.g., cardiopulmonary resuscitation); and 4) other (e.g., cardiac pacing) were recorded. For NICU patients, gestational age, birth weight, APGAR score at 1 and 5 minutes, high risk maternal risk factors (chorioamnionitis, placental abruption, etc.) were recorded.

Model Building

Quantitative Target Modeling System (QTMS V3- Applied Multivariate Algorithms, Manorville, NY) are proprietary SAS macros used for transforming raw information into optimized "bins" or intervals to be used with SAS EMS.2 neural network model building suite. Raw physiologic information was transformed creating composite variables scanned and summarized by QTMS.

This calculated a variety of statistics for each interval found within each variable including lift from base, t-tests, chi-square analysis, indices, missing data, and confidence intervals. Significant intervals from each variable were selected depending on a specified lift from base criteria. Calculating all available model building techniques, the RULE INDUCTION approach delivered the best set of specificity and sensitivity rates within the model comparison wizard. Neonatal REALTROMINS and Rapid Response REALTROMINS were created by streaming and analyzing data from the NICU and non ICU beds, respectively, into a REALTROMINS algorithm.

10

Results

Twenty-eight pediatric patients matched on age, sex, and diagnosis, half of whom died, were used for the analytic set for the neural network building process.

PICU Admission Data: The first four hours of data from the 28 matched patients were partitioned out for analyses and outcome (DIED vs. SURVIVE) recorded. This resulted in an analytic data set of 3360 records (1680 SURVIVE packets and 1680 DIED packets). Each patient thus produced 120 records for the first four hours at two minute intervals. The training data set utilized 80% of the data packets with a random 20% hold out sample of packets used to validate. The predictive capabilities of the model are demonstrated below. Identical results were found in the validation data set.

20

	Actual Patient Outcome		Total
	Survivor	Died	
Predicted to Survive	1329	21	1350
Predicted to Die	15	1322	1337
Total	1344	1343	2687

25

Predictive Capability	
Sensitivity	98.44%
Specificity	98.88%
False Positive	1.12%
False Negative	1.56%

The PICU Length of Stay (LOS): This set was analyzed representing 57,000 REALTROMINS scores with LOS from 0.29 to 12.73 days with a cumulative

length of stay of 79.41 days. The predictive capabilities of the model are demonstrated below. Six of 14 survivors and 4/14 patient deaths had 100% correct classification during the PICU LOS.

	Actual Patient Outcome		
	Survivor	Died	Total
Predicted to Survive	17,248	5,636	22,884
Predicted to Die	6,936	27,354	34,290
Total	24,184	32,990	57,174

Predictive Capability	
Sensitivity	82.92%
Specificity	71.32%
False Positive	20.23%
False Negative	24.63%

Survivor Validation Dataset: An additional 16 survivors were then analyzed during their PICU LOS representing about 25,000 REALTROMINS Scores with LOS from 0.23 to 6.36 days with 35.4 cumulative patient days. The predictive capabilities of the model were nearly identical to the PICU LOS data above.

Figure 2 is a block diagram illustrating an exemplary system for evaluating a hospital patient's risk of mortality according to an embodiment of the subject matter described herein. Referring to Figure 2, the system includes a data collection module **200** that receives inputs, referred to as feeds in Figure 2. These inputs may be physiologic signals generated by patient monitors, physiologic signals of organ function, and demographic information for a patient. Data collection module **200** may separate the data into time stamped packets. For example, where the data comprises signals output from an electrocardiogram, each data packet may include a time stamp and a peak, representing a peak ECG voltage. Messaging module **202** communicates the data packets to a real-time patient monitoring system **204**. Real-time patient monitoring system **204** includes transactional store **206** where high and low resolution data packets are stored, patient integration module **208**, which segregates the data packets into groups for individual patients and a signal processing module **210**, which generates measures of variability of at least

some of the physiologic signals. As will be described in more detail below, signal processing module **210** may receive ECG signals and use a covariance function to remove corruption from the ECG signals.

Real-time patient monitoring system **204** may also include a scoring module **212** which generates a score that indicates a patient's risk of mortality or likelihood of survival. In one exemplary implementation, scoring module **212** may determine whether a particular physiologic or demographic variable for a patient falls within a critical interval for variable that indicates that the value is predictive of mortality or likelihood of survival. Examples of critical intervals will be described below. Each time a value of a physiologic or demographic variable for the patient falls within a critical interval, scoring module **212** may record the occurrence of an event for the patient. Scoring module **212** may count the number of events for a patient over a time period, such as over a day. An output module **214** may generate, based on the count, output perceptible to a human user that indicates the patient's risk of mortality or likelihood of survival. For example, output module **214** may generate a graph that indicates a patient's score during a particular time interval. Output module **214** may transmit the graph or other indication to a physician. In one implementation, output module **214** may email a copy of the graph for a given patient to the patient's physician.

Also illustrated in Figure 2 are model development module **216** and model validation module **218**. Modules **216** and **218** may continuously update the variables used in to generate the patient's scores and may update the statistical models associated therewith to improve predictive accuracy. A data archiving module **220** may archive data for patients for further model development or validation. A system monitor **222** monitors the overall activity of the system illustrated in Figure 2.

Figure 3 is a flow chart illustrating exemplary overall steps for evaluating a hospital patient's risk of mortality according to an embodiment of the subject matter described herein. Referring to Figure 3, in step **300**, data from physiologic signals generated by patient monitors, physiologic signals of organ function, and demographic information are collected for a patient. For example, as illustrated in Figure 2, data collection module **200** may receive

feeds from patient monitors and may also receive lab data and patient demographic information. In step **302**, a measure of variability of at least one of the physiologic signals is determined. In Figure 2, signal processing module **210** may compute derivatives of a signal's output from patient monitors as one
5 measure of variability.

In step **304**, it was determined whether a value for a particular physiologic or demographic variable falls within a critical interval that indicates that the value is predictive of mortality or likelihood of surviving. For example, scoring module **212** illustrated in Figure 2 may determine whether a measured
10 or calculated value for a variable for a patient falls within a range, referred to as a critical interval, where the range indicates that the patient is likely to die or to survive. In step **306**, each time the value of for a physiologic or demographic variable falls within a critical interval, the occurrence of an event is recorded for the patient. For example, scoring module **212** may record the occurrence of
15 events of an individual patient or an individual organ in a patient. In step **308**, the number of events for a patient that occur within a time period is counted. For example, scoring module **212** may count the events for a particular patient over a time period, such as an hour or a day. In step **310**, output perceptible by a human user is generated based on the counts. The output indicates the
20 patient's risk of mortality or likelihood of surviving. For example, output module **214** illustrated in Figure 2 may generate a graph or a table that indicates the number of events for a patient. This information may be delivered to the patient's physician so that the physician can determine whether intervention is needed.

25 As set forth above, one implementation of the subject matter described herein is referred to a REALTROMINS. REALTROMINS is a medical device that analyzes in real time 1) continuous ECG signals, 2) commonly analyzed serum lab tests, and 3) diagnosis and demographic related variables. These variables or input into advanced statistical algorithms where the values are
30 compared to critical intervals to determine whether they indicate a patient is likely to die or to survive.

Time dimensioned transformations using this information included aggregating all raw data into a common time frequency as HRV information

exist in both sub-seconds and within a set of 128 beats while lab tests are intermittent based upon physician requests. First and second order derivatives, centered moving averages and standard deviations are calculated with the goal of data mining to construct a mathematical algorithm that captures viable representations of existing phenomena hidden within this database. Using QTMS V4 - Optimized Binning, a series of SAS Macros, critical intervals found with the total distribution of each variable were calculated. These bins are selected base upon their ability to lift the mortality rate from base. This is done to focus attention onto those specific ranges with significantly increased mortality rates and represent them as a binary inclusionary condition. These viable representations must be sufficiently robust, based upon their "parameter-estimation", so as to repeat any predicted classification onto to an independent hold out sample. Different classes of algorithms were found within SAS-EM5.3 (Statistical Analysis Systems, Cary NC USA) ranging from standard regression, to rule induction, and several neural networks.

One aim of the subject matter described is to provide REALTROMINS bioinformatics systems for use in hospitalized children (PICU) in the ICU and Peds Rapid Response (RR).

Classes of variables in REALTROMINS.

The following classes of variables may be used by the system illustrated in Figure 2 to evaluate a patient's risk of mortality.

- 1) Continuous physiologic signals (Le. electrocardiogram (ECG));
- 2) Advanced measures of variability of the physiologic signal (i.e. power spectral analysis);
- 3) Physiologic based measures of organ function important in predicting mortality (i.e. serum glucose) and measures of their variability
- 4) Demographic (i.e. age) and diagnosis related (i.e. post-operative status) predictors of mortality

Specific Data Elements.

The following physiologic based measures of organ function may be used by the system illustrated in Figure 2 to evaluate a hospital patient's risk of mortality.

- 5 **Fluids:** sodium, potassium, glucose, urea nitrogen, creatinine, ionized calcium, magnesium, phosphorus, whole blood glucose measurements, urine output, fluids administered (IV and/or enteral), weight
 Hematologic: hemoglobin, platelet count, prothrombin time-INR, partial thromboplastin time, D-dimers;
- 10 **Hepatic:** bilirubin, transaminases, albumin, total protein, amylase, lipase;
 Cardiac: heart rate, noninvasive blood pressure (systolic and diastolic)*, troponin, CK- MB, lactate level;
 Respiratory: rate, pH, PCO₂, PO₂, HC0₃, O₂ saturation via pulse oximetry, fractional inspired oxygen
- 15 **Neurologic:** Glasgow coma score, pain score;
 Immunologic: temperature, white cell count, absolute neutrophil and lymphocyte count, c reactive protein.

Demographic and Diagnosis Related Predictors:

- 20 The following demographic variables may be used by the system illustrated in Figure 2 to evaluate a hospital patient's risk of mortality.
 1) Age; 2) General diagnostic groups (i.e. oncologic disease); and 4) Other (nursing acuity score, therapeutic intervention score, etc.

- 25 One previous version of REALTROMINS that did not use critical intervals has been in operation since July 2008 and has expanded to include all 130 beds in the NC Children's Hospital. In addition, all of the low resolution data streams from the entire 729 bed adult hospital (excluding the monitor waveform and numeric data) are now accessible to our research environment. Since July 2008, over 2,400 pediatric patient records have been fully acquired,
- 30 processed, and archived. Uniquely, the system is created and optimized for real-time clinical decision support research, unlike conventional clinical information systems, that are either geared towards clinical operations or retrospective data collection and review. This feature is critical for

understanding and overcoming production constraints on clinical decision support research.

Continuous ECG Data: Streaming ECG Data was captured at 224 data points per second tagged with a UTC millisecond time stamp provided by the ICU monitoring system (SpaceLabs Medical). To create time-domain variables, the R peaks (fiducial point) had to be accurately identified. While the monitoring system has an analog process for estimating the location of R peaks, this process does not always identify the correct number or location of the R peak. Figures 4A and 4B illustrate the raw ECG signal with peaks marked using SpaceLabs algorithms.

Stripping out the actual Rpeak time stamp and then counting binned-values of the time difference between peaks in increments of 100 milliseconds calculated time domain variables. While current research had focused only on the pNN50 value we found significance elsewhere among the range of this distribution. Significance was found among different lags, e.g., pNN700, pNN800, and pNN400. Figures 5A and 5B are tables illustrating means and standard deviations for pNNXXX distributions (dead1=arrested, 0=alive). From the NN interval data (time between successive R waves), SDNN (standard deviation of the NN (consecutive normal sinus intervals), RMSSD (root mean square successive difference), NN50 (the number of NN intervals with differences > 50 ms), pNN50 (the proportion of NN50 intervals divided by the total NN intervals), pNNx series (where x=time>0 ms), NN range, coefficient of variation, and the min/max heart period for each 128 beat minor time epoch are calculated.

Spectral data was extracted by stripping on the fixed 64 records for each group of records. Frequency bands were grouped as follows using the interpolated heart rate spectra values to calculate area within that specific band of frequencies. These frequencies corresponded to the ultra-low frequency (ULF), very low frequency (VLF), low frequency (LF), and high frequency (HF) bands.

Transformations were performed upon the interpolated heart rate spectra values to convert them into variance. The variance present in the HR data is obtained from the power spectrum according to methods described in

Bendat JS, Piersol AG. Engineering Applications of Correlation and Spectral Analysis. 2nd Edition. New York: John Wiley & Sons, Inc.; 1993, the disclosure of which is incorporated herein by reference in its entirety. Since our HR data was normalized, the mean is zero, and the sum of the power in the frequency

5 bands is directly equal to the variance. We chose to generate our spectra using the HR time series data in a similar manner to methods used by Pomeranz B, MacAulay R, Caudill M. Assessment of autonomic function in humans by heart rate spectral analysis. *Am J Physiol* 1985;248:H151-3, the disclosure of which is incorporated herein by reference in its entirety. Because

10 we use the normalized HR data and not the heart beat period time series, our units of power are *bpm2/bpm2* (unitless) as opposed to *ms2* or normalized units reported by some other investigators. The value for power reported within each frequency band will vary, depending on the units, but the ratios of interest (e.g. low frequency power to high frequency power, LF/HF) will not be affected.

15 Figures 6A and 6B are tables illustrating means and standard deviations from spectral distributions where dead1=arrested, 0=alive.

It can be seen from Figures 6A and 6B that the standard deviations exceed the means indicating nonlinear distributions. Using QTMS optimal bins, critical intervals were developed by isolating specific sub-ranges within the total

20 distribution of each variable. They are color coded by importance in identifying dead patients by focusing attention to cut points masked by the means. Standard means testing focuses on values of .54 for the dead while the extreme bin below .34 indexes at 344. This approach is used on every variable to recode new variables with significant lift. Figure 7 is a table illustrating

25 critical intervals used by the system illustrated in Figure 2 to evaluate a patient's risk of mortality. In Figure 7, intervals 1 and 2 may be classified as critical intervals because they have mortality rates of 47.63 and 20.53%, respectively. While not accurately predictive of mortality alone, a variable that falls in intervals 1 and 2 can be considered in a model with other variables that

30 can be highly predictive of a patient's risk of mortality. Exemplary models that combine variables will be set forth in detail below.

Eight additional transformations to each raw variable were also created and include centered moving averages, centered moving standard deviations

nested across first and second order derivatives in attempts to provide information on how these values change over the course of each patient's stay. Figures 8A and 8B illustrate transformable variables that may be used by the system illustrated in Figure 2 to predict a hospital patient's risk of mortality.

5 Physiologic based measures of organ function important in predicting mortality: were defined by a battery of 96 laboratory and streaming physiological tests collected during the patients stay in the PICU. These tests were collected based upon clinical need which varied by patient and varied over time within patient. This was accomplished by looking within each
10 variable's distributions and identifying intervals that index high against base-mortality using the optimized binning procedure found within QTMS V4.0 and previously described. This procedure was also performed on the eight additional transformations for each variable (centered moving averages, standard deviation on first and second derivatives). Figures 9A and 9B
15 illustrate tables that contain laboratory tests that are available for evaluating a hospital patient's risk of mortality.

 The basic idea is to find intervals within the total distribution where the target distribution (arrested patients) occurs in different percentages. The output found in below helps identify these intervals by first binning the entire
20 distribution into deciles or 10 equaled sized partitions. These intervals are used to calculate a mortality rate within that range and index it to the overall mortality rate. The table below that displays those intervals where the index is greater than 120. From the 98 raw labs and streaming physiologic data, over 500 composite and transformed variables were reported on with 2,400 intervals
25 from that set indexing 120+. Figures 10A and 10B are tables illustrating composite and transformation variables that may be used by the system illustrated in Figure 2 to evaluate a hospital patient's risk of mortality. Figures 11A and 11B are tables illustrating critical intervals for the variables of Figures 10A and 10B that may be used to evaluate a hospital patient's risk of mortality.

30 The advantage to using these quantitatively derived intervals is that any means testing tends to focus attention on incorrect areas within the total distributions when the data is not normally distributed. A t-test for Sodium compares a mean of 141 for live patients versus 137 for arresting patients.

You find flat indexes about those values in the graphic below. The data is not normal but skewed. By recoding the green intervals into a (0,1) we intensify the strength of the effect. We also have an advantage in that we can now include patients who did not have any test by coding them as a 0. Figures 12A and 12B are respectively a table and a graph of critical intervals for sodium content of live and dead patients. Figures 13A and 13B are graphs of sodium frequency versus time. Figures 14A - 14C illustrate event counts for different organs for different patients. Figures 15A - 15D illustrate critical intervals that may be used by the system illustrated in Figure 2 for evaluating a patient's risk of mortality. The critical intervals in Figures 15A - 15D may be determined using a single factor scan procedure found with QTMSV3.2 profiling for extreme derogatory ranges within the battery of lab tests.

Developing a Common Time Hierarchy. This study was an accumulation of patient data collected continuously during the PICU hospitalization. The units of time for each grouping of variables were different. The units for the time domain variables were sub-second differences between R peaks, while the time span represented in the frequency domain was based on 128 bpm). Finally the units for intermittent variables were logged at the actual clock time taken. To join these 4 groups of data, it was determined to put this information into a common time frame of 2 minutes. This was accomplished using SAS procedure Expand (part of the ETS Suite within SAS) which can combine time series with different frequencies using various interpolative methods that can be used to convert raw data into a higher frequency series or aggregate down to a lower frequency series. Figures 16A and 16B illustrate raw and extrapolated data generated using the procedure Expand. Figure 16A illustrates intermittent original data and Figure 16B illustrates extrapolated data into increments.

Following these procedures, we extrapolated all data from all sources, calibrated into 2 minute epochs for all patients. This resulted in an analytic data set of 1,000,000 two minute records that contain the 500+ variables and transformations described before. Each packet represents two minutes of clock time with information from the pNN category, spectral information, and derogatory markers from the lab tests and demographics collected.

In conclusion, we have taken data from disparate sources and collected along difference time increments and built an integrated single source analytic file all calibrated into two minute increments.

To answer issues about overfitting and the fact that many of these critical intervals are correlated, an added step was performed with the NICU model which was to reduce the number of intervals by various rank reduction techniques (PCA) and model using the factor scores. Because each factor score is continuous, a Two-Group Multiple Discriminant was used and is reported.

Multivariate Model Building: The composite file described above was analyzed using SAS Enterprise Miner 5.3 with the logistic regression model yielding absolute best fit among several close alternatives. All data through December 2009 yielded a file of 750,000 2 minute packets each containing the 500+ raw and transformed variables along with the 2,400 optimized bins (based upon the critical intervals among these 500+ variables). A 40% training sample (304,929 packets) was used to model build with a 30% validation sample (228,696 packets) and a 30% test sample (228,698 packets) held out. The summary of all samples is provided below.

Training Sensitivity = $38,739 / (38,739 + 3526) = 91.65\%$
 Training Specificity = $261,586 / (261,586 + 1078) = 99.58\%$
 Validation Sensitivity = $29,063 / (29,063 + 2636) = 91.68\%$
 Validation Specificity = $196,107 / (196,107 + 890) = 99.54\%$

The testing sample produced identical results as the validation sample concluding the model is robust.

Conclusions indicate correctly NOT alarming LIVE patients 99% of the time and correctly alarming DEAD patients 91% of the time during the entire course of their PICU stay. Conversely the model misses a DEAD packet 9% of the time. While this represents each 2 minute result, accumulating signals over the course of 30 minutes to 1 hour would yield a clearer picture and a better summarization of the patient's current status.

Interactive Reporting: As we collect new patient information and process it through the same algorithms just described, using COGNOS data management software available from IBM corporation and this test model, REALTROMINS risk of mortality score with a drill down of the individual components contributing to that score will be provided in real time allowing for diagnostic updates every 2 minutes. Figure 17A is a graph and Figure 17B is a table illustrating exemplary scores that may be generated for a patient. The scores may be presented on a per-organ and/or per-patient basis. Such scores may be emailed to a physician in graphical format or in tabular format to allow the physician to make a quick evaluation of the patient's risk of mortality.

Information visualization is to present significant and non-obvious patterns within these vast high-dimensional datasets and provide users with intuitive interface in their exploratory analysis. Many salient patterns may be only visible in subspaces of the data. Dimensionality reduction and clustering have become an effective approach to capture and organize these patterns, aid in the removal of irrelevant or redundant information, and enhances the comprehensibility of the whole dataset.

In this project, each patient is monitored continuously and a high-dimensional feature vector is derived for every 2-minute long interval. A method to organize and summarize these feature vectors so that a user can navigate through them easily is desirable. In order to generate a consistent yet succinct view of all data, we need to establish and evaluate the relationships between these feature vectors. Given the pairwise relationships between feature vectors, a representation that can efficiently organize them is needed. This includes the ability to group and summarize similar vectors and the ability to allow users to examine the relationship between feature vectors and/or groups of feature vectors at various resolutions seamlessly. Ideally, we want to derive plausible hierarchical categorizations of feature vectors from their pairwise dissimilarities. We plan to project these high dimensional vectors onto some low dimensional space and apply hierarchical clustering.

We have continued to improve upon our existing reporting technologies imbedded in REALTROMINS. Exemplary unconventional and exceptionally innovative aspects to our approach are listed below:

- streaming real time continuous physiologic data for the entire length of ICU hospitalization;
- interfacing multiple disparate forms of clinically important patient data into a common database for research development of an advanced bioinformatics system;
- providing an informatics system storing high and low resolution data to readily allow for frequent recalibrations of the predictive models;
- providing mortality predictions based on statistical principles, but utilizing novel critical interval based techniques
- providing a real time continuously updated risk of mortality assessment across a broad spectrum of patients; and
- providing a graphical interface that allows bedside care givers to understand and act upon large sets of clinically important data.

As set forth above, one aspect of the subject matter described herein uses a model that includes critical intervals for multiple variables to which monitored values for a patient are compared in order to generate an indication of the patient's risk of mortality. The following is an example of a NICU model that may be implemented by the system illustrated in Figure 2 to evaluate a patient's risk of mortality:

Variable	0	1
Constant	-0.20091	-6.81 060
Factor1	-0.28830	1.65814
Factor2	0.00604	-0.00391
Factor3	-0.01 189	0.07554
Factor4	-0.69595	4.07369
Factor5	-0.28752	1.66290
Factor6	-0.26048	1.52776
Factor7	-0.08562	0.47352
Factor8	-0.20929	1.24321
Factor9	-0.21 608	1.251 33
Factor1 0	-0. 13746	0.80579
Factor1 1	-0. 15405	0.90698
Factor1 2	-0.06925	0.38799
Factor1 3	-0.0991 3	0.57833
Factor1 4	-0.0871 1	0.531 04
Factor1 5	-0.08240	0.53599
Factor1 6	-0. 12089	0.71 464
Factor1 7	-0.09577	0.54900
Factor1 8	-0. 141 58	0.821 08
Factor1 9	-0. 13160	0.74267
Factor20	-0. 1361 3	0.80762
Factor21	-0.08358	0.51459
Factor22	-0.06739	0.41 867
Factor23	-0. 10 166	0.60537

	Factor24	-0.18174	1.05345
	Factor25	-0.23663	1.36195
	Factor26	-0.01213	0.10820
	Factor27	0.05107	-0.29622
5	Factor28	-0.14207	0.78863
	Factor29	-0.06722	0.40433
	Factor30	-0.06454	0.38491
	Factor31	-0.07605	0.42829
	Factor32	-0.11187	0.67386
10	Factor33	-0.03017	0.13137
	Factor34	-0.01941	0.13086
	Factor35	-0.02694	0.14742
	Factor36	-0.03483	0.18450
	Factor37	-0.04668	0.29202
15	Factor38	-0.03601	0.20023
	Factor39	0.05874	-0.31724
	Factor40	-0.06479	0.39679
	Factor41	-0.09578	0.55065
	Factor42	-0.07738	0.45175
20	Factor43	-0.05947	0.30176
	Factor44	-0.16494	0.96125
	Factor45	-0.04428	0.24945
	Factor46	-0.02387	0.15645
	Factor47	-0.05996	0.33865
25	Factor48	-0.10723	0.61003
	Factor49	-0.03515	0.18770
	Factor50	-0.02083	0.15172
	Factor51	-0.04029	0.22747
	Factor52	-0.06045	0.36015
30	Factor53	-0.04780	0.30002
	Factor54	-0.07062	0.39624
	Factor55	-0.07177	0.41246
	Factor56	-0.05444	0.29256
	Factor57	-0.05473	0.35851
35	Factor58	-0.05819	0.36887
	Factor59	-0.05257	0.29494
	Factor60	-0.04006	0.25405
	Factor61	-0.05685	0.28521
	Factor62	-0.05115	0.29318
40	Factor63	-0.05635	0.34883
	Factor64	-0.06668	0.39700
	Factor65	-0.06095	0.36169
	Factor66	-0.05292	0.28978
	Factor67	-0.04499	0.21666
45	Factor68	-0.00444	0.01509
	Factor69	-0.04559	0.23776
	Factor70	-0.04526	0.26049
	Factor71	-0.11500	0.65056
	Factor72	-0.03565	0.19519
50	Factor73	-0.04075	0.22143
	Factor74	-0.03697	0.21035
	Factor75	-0.03543	0.21661
	Factor76	-0.04793	0.27253
	Factor77	-0.05694	0.29992
55	Factor78	-0.08618	0.50202
	Factor79	-0.03524	0.20943
	Factor80	-0.04470	0.24862
	Factor81	-0.05235	0.29493
	Factor82	-0.07565	0.44194
60	Factor83	-0.03100	0.18079

	Factor84	-0.04348	0.25106
	Factor85	-0.05120	0.28331
	Factor86	-0.03110	0.18381
	Factor87	-0.03059	0.15312
5	Factor88	-0.04195	0.21911
	Factor89	-0.02403	0.13863
	Factor90	-0.03540	0.18019
	Factor91	-0.08765	0.51857
	Factor92	-0.04127	0.22865
10	Factor93	-0.01852	0.13274
	Factor94	-0.04286	0.24996
	Factor95	-0.00759	0.01358
	Factor96	-0.06443	0.35972
	Factor97	-0.03615	0.15639
15	Factor98	-0.04045	0.25249
	Factor99	0.00357	-0.03108

In the model, each factor represents a physiologic variable or combination of such variables. The "0" column represents beta weights for Live NICU patients while the "1" column represents beta weights for the dead NICU patients. In addition we have also produced PMML formatted scoring algorithms of the PICU and RR population. Abbreviated examples can be found for PICU in Figure 18A and for RR in Figure 18B. Although different critical intervals were identified for each patient population (NICU, PICU and RR), all follow the same reporting outlined in Figures 17A- 17C. Complete versions of these models can be found in the above-referenced provisional patent application.

Improving ECG Signal Integrity: Pre-processing of the ECG:

According to another aspect of the subject matter described herein, the predictive accuracy of the models can be improved by pre-processing the ECG signal using a covariance function to lessen the effects of corruption in the ECG signal. Traditional methods of "R" wave detection of the ECG use a combination of signal conditioning, QRS detection and post-processing. QRS detection methods include amplitude threshold, first derivative threshold, moving average filters, high pass filters, cross correlation and others. Our approach is to use the covariance function (formula 1) in conjunction with statistical filtering to detect the "R" wave. We have found the covariance function is superior for QRS detection than standard cross-correlation or correlation coefficient function, since it removes baseline shifts and does not amplify other portions of the ECG other than the "R" wave. Figure 19 is an ECG

signal which has been corrupted: from second 2-4 with baseline shift due to motion artifact from deep breathing; from seconds 3-7 with moderate EMG artifacts due to a Valsalva maneuver; and from seconds 7-11 with intense pectoralis major muscle flexion.

5

$$C_{xy}(m) = \frac{1}{N} \sum_{n=0}^N (x(n) - \mu_x)(y(n+m) - \mu_y) \quad (\text{formula 1})$$

Using the covariance function (formula 1), the corruption can be removed. For example, when missing data is seen in ECG signals we "fill-in" using the correlation (covariance) with secondary sources (e.g., respiration) to estimate the missing ECG value. When high covariance(correlation) with secondary feeds exist, this approach has yielded significant value over a randomly assigned/blank values. In addition, when looking across multiple secondary feeds simultaneously we have reported more accurate estimates (fill-in value) then when just using single source.

The covariance between the QRS template (Figure 20) and the contaminated ECG (Figure 19) is determined; with the output shown in Figure 21.

20 Post Processing

Other investigators have used finite impulse response (FIR), infinite impulse response (IIR) filters, adaptive filtering, cubic spline fitting, moving average filtering, median filtering, weighted average filtering, mathematical morphology and neural networks to aid in the identification of the QRS wave. We have implemented a quadrature-mirrored 2nd order linear high-pass IIR filter to remove the baseline shifts. The output is shown in Figure 22.

The R wave was then distinguished using a 2-sigma statistical filter, followed by a low-pass discriminatory filter, as shown in Figure 23.

The output of the 2-sigma statistical filter is then applied to a Threshold and First Derivative filter, with the output of this final stage shown in Figure 24. In this figure, the original contaminated ECG is displayed, with the output of the

several filtering steps described above shown as binary values (1 or 0) along the bottom axis. The spikes (value 1) mark the occurrence of the "R" wave. The ECG signal in Figure 24 represents a section of a worst-case scenario test signal with a SNR of 0.5 (-6.1 dB), which only occurs occasionally with poor
5 sensor connectivity or during extreme patient activity (0.4 Hz - respiration, 3 Hz cardioballistic, 60 Hz - power line, 0-155Hz, white noise). Applying this detection procedure, we are able to correctly detect the "R" wave on this ECG test signal at 92% sensitivity, with 20% false positives.

We have tested this "R" wave detection algorithm in a research study
10 with 24 hr monitoring of sleep patients (20 subjects). While our current method is sufficient for ECG signal processing when off-line correction is appropriate, it is not sufficient for real-time continuous "R" wave detection for alarming, which must have very high real-time fidelity. High level Alarming for this condition must be prevented to achieve caregiver alertness through confidence in the
15 alarm situation. We will improve our detection methodology to 99% detection with 5% false positives under this extreme condition of ECG contamination.

Alarming

Even following signal filtering and enhanced identification of "R" waves
20 using cross-correlation of the ECG with QRS templates, waveform sections with missed beats or falsely marked "R" waves are inevitable, resulting in false alarms. We will employ a multi-tiered method for minimizing false alarms.

Heart Rate Alarms:

25 Three parallel records of heart period time-series will be derived following the initial signal conditioning and RR interval detection from the two ECG signals and the invasive or non-invasive pulse waveform obtained from the pulsed oximeter. Missing or false detections of the R-wave will be implied from unusually long (multiple of a single interval) or short (fractional) intervals.
30 Paired comparisons will be made between these three data sets and improperly detected events will be corrected when they coincide with the expected interval. When a long segment of any of the three waveforms is false, a low-level alarm will be issued (poor sensor contact, high level noise, or

sensor disconnect). When long segments (5 seconds) of all three waveforms are false, a high level alarm will be issued.

- It will be understood that various details of the subject matter described herein may be changed without departing from the scope of the subject matter
- 5 described herein. Furthermore, the foregoing description is for the purpose of illustration only, and not for the purpose of limitation.

CLAIMS

What is claimed is:

1. A method for evaluating a hospital patient's risk of mortality, the
5 method comprising:
 - collecting data from physiologic signals generated by patient
monitors, physiologic signals of organ function, and demographic
information for a patient;
 - 10 determining a measure of the variability of at least one of the
physiologic signals;
 - analyzing the data and the measure of variability to determine
whether value for a particular physiologic or demographic variable falls
within a critical interval for the variable that indicates that the value is
predictive of mortality or survival;
 - 15 each time a value for a physiologic or demographic variable for
the patient falls within a critical interval, recording the occurrence of an
event for the patient;
 - counting the number of events for the patient over a time period;
and
 - 20 generating, based on the count, output perceptible by human
user that indicates the patient's risk of mortality or likelihood of survival.
2. The method of claim 1 wherein collecting the data includes
separating the data into time stamped packets.
3. The method of claim 2 wherein the data comprises data from an
25 electronic organ monitor, wherein separating the data into time
stamped packets includes identifying a peak value and a time stamp
for each packet.
4. The method of claim 3 wherein determining a measure of the
variability comprises computing derivatives in the data from the peak
30 values and the time stamps.
5. The method of claim 1 wherein each critical interval comprises a
range of values for which mortality or non-mortality rate for known
patient outcomes is above a threshold.

6. The method of claim 1 wherein counting the number of events for the patient includes counting the number of total events for the patient in a day.
7. The method of claim 1 wherein counting the number of events includes counting the number of events per organ for the patient.
8. The method of claim 1 wherein generating the output includes generating a graph of event counts for the patient.
9. The method of claim 1 comprising electronically transmitting the graph to a physician.
10. The method of claim 1 wherein the patient comprises one of a neonatal intensive care unit patient, a pediatric intensive care unit patient, or a patient in the hospital prior to admission to an intensive care unit.
11. The method of claim 1 wherein the physiologic signals include an electrocardiogram (ECG) signal and wherein the method further comprises using a covariance function to reduce corruption in the ECG signal.
12. A system for evaluating a hospital patient's risk of mortality, the system comprising:
 - a data collection module for collecting data from physiologic signals generated by patient monitors, physiologic measures of organ function, and demographic information for a patient in a hospital;
 - a signal processing module for computing a measure of the variability of the physiologic signals;
 - a scoring module for analyzing the data and the measure of variability to determine whether a value for a particular physiological or demographic variable falls within a critical interval for the variable that indicates that the value is predictive of mortality or survival, for, each time the value for a physiologic or demographic variable for the patient falls within a critical interval, recording the occurrence of an event for the patient and counting the number of events for the patient over a time period; and

an output module for generating, based on the count, output perceptible by human user that indicates the patient's risk of mortality or likelihood of survival.

- 5 13. The system of claim 12 wherein the data collection module is configured to separate the data into time stamped packets.
14. The system of claim 13 wherein the data comprises data from an electronic organ monitor, wherein separating the data into time stamped packets includes identifying a peak value and a time stamp for each packet.
- 10 15. The system of claim 14 wherein the signal processing module is configured to compute derivatives in the data from the peak values and the time stamps.
16. The system of claim 12 wherein each critical interval comprises a range of values for which mortality or survival rate for known patient outcomes is above a threshold.
- 15 17. The system of claim 12 wherein the scoring module is configured to count the number of total events for the patient in a day.
18. The system of claim 12 wherein the scoring module is configured to count the number of events per organ for the patient.
- 20 19. The system of claim 12 wherein the output module is configured to generate a graph of event counts for the patient.
20. The system of claim 12 wherein the output module is configured to transmit the graph to a physician.
21. The system of claim 12 wherein the patient comprises one of a neonatal intensive care unit patient, a pediatric intensive care unit patient, or a patient in the hospital prior to admission to an intensive care unit.
- 25 22. The system of claim 12 wherein the signal processing module is configured to receive and process an electrocardiogram (ECG) signal and to use a covariance function to reduce corruption in the ECG signal.
- 30

23. A non-transitory computer readable medium having stored thereon executable instructions that when executed by the processor of a computer control the computer to perform steps comprising:

5 collecting data from physiologic signals generated by patient monitors, physiologic signals of organ function, and demographic information for a patient;

 determining a measure of the variability of at least one of the physiologic signals;


10 analyzing the data and the measure of variability to determine whether value for a particular physiologic or demographic variable falls within a critical interval for the variable that indicates that the value is predictive of mortality or survival;

15 each time a value for a physiologic or demographic variable for the patient falls within a critical interval, recording the occurrence of an event for the patient;

 counting the number of events for the patient over a time period; and

 generating, based on the count, output perceptible by human user that indicates the patient's risk of mortality or likelihood of survival.

20







FIG. 2

3/32

**FIG. 3**




FIG. 4A

5/32




FIG. 4B

6/32

**DEAD PICU
THE MEANS PROCEDURE**

VARIABLE	MEAN	Std Dev	VARIANCE	MINIMUM	MAXIMUM	N
pnn50	0	0	0	0	0	105664
pnn100	0	0	0	0	0	105664
pnn200	0.0071923	0.0528221	0.0027902	0	0.9922780	105664
pnn300	0.1080836	0.2913443	0.0848815	0	0.9937888	105664
pnn400	0.3028720	0.4268584	0.1822081	0	0.9935897	105664
pnn500	0.3405127	0.4214867	0.1776511	0	0.9934211	105664
pnn600	0.1618855	0.3109868	0.0967128	0	0.9934211	105664
pnn700	0.0234902	0.1100181	0.0121040	0	0.9923077	105664
pnn800	0.0153742	0.0767250	0.0058867	0	0.9923077	105664
pnn900	0.0084066	0.0456894	0.0020875	0	0.9534884	105664
pnn1000	0.0057024	0.0335443	0.0011252	0	0.7224490	105664
pnn1000o	0.0173942	0.0710719	0.0050512	0	0.9840000	105664

FIG. 5A

**LIVE PICU
THE MEANS PROCEDURE**

VARIABLE	MEAN	Std Dev	VARIANCE	MINIMUM	MAXIMUM	N
pnn50	6.1405662E-9	4.9759771E-6	2.476035E-11	0	0.0040323	656659
pnn100	0.000018966	0.0036391	0.000013243	0	0.7142857	656659
pnn200	0.000500129	0.0192136	0.000363164	0	0.9922481	656659
pnn300	0.0049646	0.0600159	0.0036019	0	0.9939394	656659
pnn400	0.1388401	0.3025583	0.0915415	0	0.9939759	656659
pnn500	0.3566494	0.4133237	0.1708365	0	0.9934211	656659
pnn600	0.2273262	0.3440774	0.1183802	0	0.9930556	656659
pnn700	0.1223216	0.2496523	0.0623263	0	0.9926471	656659
pnn800	0.0629374	0.1584039	0.0250918	0	0.9923077	656659
pnn900	0.0326490	0.0938397	0.0083059	0	0.9922481	656659
pnn1000	0.0135719	0.0488648	0.0023875	0	0.9397590	656659
pnn1000o	0.0235014	0.0863901	0.0074632	0	0.9920000	656659

FIG. 5B

7/32

**DEAD PICU
THE MEANS PROCEDURE**

VARIABLE	MEAN	Std Dev	VARIANCE	MINIMUM	MAXIMUM	N
RATIO	0.5475689	0.3979583	0.1583708	0.038972	10.8972262	105664
ULF	9.7039194E-6	0.000281968	7.950572E-8	1.715726E-10	0.0540024	105664
VLF	0.000890117	0.0329209	0.0010838	1.3612728E-6	6.8905225	105664
LF	0.0020710	0.0785836	0.0061754	3.1578073E-6	17.5578519	105664
HF	0.0032880	0.0202568	0.000410338	0	4.3284978	105664
VHF	0.0054152	0.0114744	0.000131661	0	1.5805266	105664

FIG. 6A

**LIVE PICU
THE MEANS PROCEDURE**

VARIABLE	MEAN	Std Dev	VARIANCE	MINIMUM	MAXIMUM	N
RATIO	0.8650416	3.3508340	11.2280887	0.0336705	705.7305467	656659
ULF	0.000056382	0.0056345	0.000031748	7.6014E-10	4.3003126	656659
VLF	0.0035131	0.2877154	0.0827801	1.7601426E-6	190.5933976	656659
LF	0.0044045	0.1710344	0.0292528	0.000014352	58.6030532	656659
HF	0.0044401	0.0888781	0.0078638	0	50.0547221	656659
VHF	0.0047622	0.0447621	0.0020036	0	28.7881662	656659

FIG. 6B

8/32

ENTER VAL	RANGE	PICU COUNT	DEAD COUNT	INDEX	Prob	LEFT Conf	MORTALITY RATE	RIGHT Conf	SELECT
1	0.03367047-0.34640151	76,241	36,311	344	0.000	47.27	47.63	47.98	HIGH
2	0.346402112-0.439258343	76,232	15,652	148	0.000	20.25	20.53	20.82	HIGH
3	0.439258833-0.49648596	76,232	12,479	118	0.000	16.11	16.37	16.63	FLAT
4	0.49648620-0.544375899	76,232	8,037	85	0.000	11.50	11.72	11.95	FLAT
5	0.544377051-0.612521286	76,232	6,506	62	0.000	8.34	8.53	8.73	LOW
6	0.61252217-0.715749856	76,232	5,249	50	0.000	6.71	6.89	7.07	LOW
7	0.715750008-0.861432643	77,465	5,275	49	0.000	6.63	6.81	6.99	LOW
8	0.861448546-1.119730555	76,233	5,878	56	0.000	7.52	7.71	7.90	LOW
9	1.119748978-1.499777383	76,233	5,397	56	0.000	7.55	7.71	7.93	LOW
10	1.499780634-795.7395467	74,991	3,480	33	0.000	4.49	44.6	4.79	LOW

FIG. 7

9/32

**DEAD PICU
THE MEANS PROCEDURE**

VARIABLE	MEAN	Std Dev	VARIANCE	MINIMUM	MAXIMUM	N
RATIO	0.5475689	0.3979583	0.1583708	0.0381972	10.8972262	105664
RATIO_CMOVAV	0.5475493	0.2778088	0.0771777	0.1291440	2.6945082	105664
RATIO_CMOVSTD	0.2124678	0.1931013	0.0372881	-1.73472E-18	3.8410779	105664
DERIV1_RATIO	9.6549427E-8	0.0023481	5.5135594E-6	-0.0551619	0.0559047	105664
DERIV1_RATIO_CMOVAV	6.6880406E-8	0.000063608	4.0459309E-9	-0.0015887	0.0015779	105664
DERIV1_RATIO_CMOVSTD	0.0018233	0.0015016	2.2547139E-6	-1.43821E-93	0.0120505	105664
DERIV2_RATIO	-5.297347E-9	0.000028850	8.3234E-10	-0.0012118	0.000904234	105664
DERIV2_RATIO_CMOVAV	-3.199729E-9	9.2077017E-7	8.478177E-13	-0.000035571	0.000028564	105664
DERIV2_RATIO_CMOVSTD	0.000022039	0.000018609	3.463014E-10	-1.43821E-93	0.000226397	105664

FIG. 8A

**LIVE PICU
THE MEANS PROCEDURE**

VARIABLE	MEAN	Std Dev	VARIANCE	MINIMUM	MAXIMUM	N
RATIO	0.8650416	3.3508340	11.2280887	0.0336705	795.7395467	656659
RATIO_CMOVAV	0.8648648	1.4305640	2.0465134	0.1772132	147.5979595	656659
RATIO_CMOVSTD	0.3853615	3.0298718	9.1801231	-4.85723E-17	310.3030093	656659
DERIV1_RATIO	3.7340002E-7	0.0114227	0.000130479	-4.1634044	4.2015571	656659
DERIV1_RATIO_CMOVAV	3.1721199E-7	0.000654915	4.289138E-7	-0.1252407	0.1252371	656659
DERIV1_RATIO_CMOVSTD	0.0028854	0.0111235	0.000123732	-3.35583E-93	1.303643	656659
DERIV2_RATIO	4.5103394E-9	0.000118886	1.4133862E-8	-0.0384450	0.0384463	656659
DERIV2_RATIO_CMOVAV	3.7799388E-9	2.6787591E-6	7.17575E-12	-0.000797139	0.000797070	656659
DERIV2_RATIO_CMOVSTD	0.000034288	0.000114582	1.3129015E-8	-3.35583E-93	0.0115419	656659

FIG. 8B

10/32

ABSOLUTE_NEUT_CT	FIBRINDGEN-F1	O2 SATURATION	PT
ALT	CGT	O2 SATURATION-VEN	PTH
APTT	GLUCOSE POC	PCO2	SODIUM
AST	GLUCOSE RANDOM	PCO2/POC	SODIUM (BG/CC)
BASE BALANCE	GLUCOSE WHL BLOOD	PCO2-VENOUS	SODIUM-BG/POC
BASE BALANCE/POC	GLUCOSE-BG/POC	PCO2-VENOUS/POC	TROPONIN T
BASE BALANCE-VEN	HEMATOCRIT	PH	UREA NITROGEN
BASE BALANCE-VENOUS/POC	HEMATOCRIT-BG/P	PH/POC	URIC ACID
BICARB/POC	HEMOGLOBIN	PHOSPHORUS	WBC
BICARBONATE	HEMOGLOBIN(BG/CC)	PH-VENOUS	
BICARBONATE-VEN	INR	PH-VENOUS/POC	
BICARB-VEN/POC	IONIZED CA VENOUS	PLATELET COUNT	
BILIRUBIN TOTAL	IONIZED CA-SG/POC	P02	
CALCIUM	IONIZED CALCIUM	P02/POC	
CHLORIDE	LACTATE ARTERIAL	P02-VENOUS	
CK-MB	LACTATE VENOUS	P02-VENOUS/POC	
CO2	LACTATE-PCC	POTASSIUM	
CREATINE KINASE	LACTATE-VEN/POC	POTASSIUM (BG/CC)	
CREATININE	MAGNESIUM	POTASSIUM-BG/POC	
D-DIMER QUANT	NEUT LEFT SHIFT	PROTEIN TOTAL	

FIG. 9A

ART DIAS	SPO2
ART MEAN	SpO2 PULSE RATE
ART PULSE RATE	TEMP 2
ART SYS	TEMPERATURE
CAPNO RESP RATE	UA DIAS
CORE TEMP	UA MEAN
CVP	UA SYS
ETCO2	UAP PULSE RATE
HEART RATE	UV DIAS
HR ECG 2	UV SYS
ICP	UVP MEAN
MAP	
NIBP DICS	
NIBP MEAN	
NIBP SYS	
RESP RATE	

FIG. 9B

11/32

QTMS 4.0		SINGLE FACTOR SCAN SUMMARY TABLE OF PROJECT TESTRUN THE VARIABLES ARE LISTED IN ALPHABETICAL ORDER						
#	VARIABLE	NO. OF INTERVALS	CHI- SQUARE	CHI-SQ. PROB.	T. TEST	T. TEST PROB.	MEAN OF NON.RESP.	MEAN OF RESPONDERS
1	ABSOLUTE_NEUT_CUT	10	101786.26	0.0000	-68.0354	0.0000	7.746	10.525
2	ABSOLUTE_NEUT_CUT_CMOVAV	10	169838.31	0.0000	-114.3353	0.0000	6.767	10.967
3	ABSOLUTE_NEUT_CUT_CMOVSTD	10	135636.67	0.0000	-87.2072	0.0000	3.128	4.488
4	ALT	10	67695.334	0.0000	-18.6174	0.0000	123.138	139.497
5	ALT_CMOVAV	10	127890.70	0.0000	-101.9348	0.0000	100.744	175.135
6	ALT_CMOVSTD	10	86256.561	0.0000	-57.7809	0.0000	65.323	98.400
7	ANDMOL	2	2.356	0.1248	4.1232	0.0000	0.000	0.000
8	APTT	10	149982.52	0.0000	-59.7890	0.0000	38.881	46.767
9	APTT_CMOVAV	10	183124.95	0.0000	-55.6134	0.0000	43.090	49.290
10	APTT_CMOVSTD	10	183231.68	0.0000	-25.1562	0.0000	12.151	13.995
11	AST	8	76281.535	0.0000	-50.7186	0.0000	221.533	457.244
12	AST_CMOVAV	8	59752.489	0.0000	-74.5702	0.0000	135.681	469.084
13	AST_CMOVSTD	8	75688.803	0.0000	-133.3939	0.0000	110.280	218.933
14	BILIRUBIN_TOTAL	8	118250.25	0.0000	-92.4685	0.0000	0.720	1.696
15	BILIRUBIN_TOTAL_CMOVAV	8	159809.90	0.0000	-210.8066	0.0000	0.662	2.390
16	BILIRUBIN_TOTAL_CMOVSTD	7	166362.96	0.0000	-252.6595	0.0000	0.205	2.029
17	CALCIUM	10	61034.470	0.0000	55.1787	0.0000	9.091	8.898
18	CALCIUM_CMOVAV	10	104352.22	0.0000	139.5541	0.0000	9.119	8.849
19	CMOVAV_CMOVSTD	10	72683.616	0.0000	22.0857	0.0000	0.644	0.615
20	CAPNO_RESP_RATE	10	45764.334	0.0000	-11.3172	0.0000	24.588	25.261
21	CAPNO_RESP_RATE_CMOVAV	10	44966.418	0.0000	-9.9354	0.0000	25.023	25.560
22	CAPNO_RESP_RATE_CMOVSTD	10	48770.078	0.0000	104.1389	0.0000	3.229	1.585
23	CHLORIDE	9	66526.244	0.0000	1.4763	0.1399	102.394	102.342
24	CHLORIDE_CMOVAV	10	93255.426	0.0000	9.8423	0.0000	102.143	101.887
25	CHLORIDE_CMOVSTD	10	65547.947	0.0000	-1.4849	0.1376	4.790	4.801
26	CO2	9	55048.530	0.0000	31.8660	0.0000	24.945	24.368
27	CO2_CMOVAV	10	138557.64	0.0000	32.0350	0.0000	24.862	24.481
28	CO2_CMOVSTD	10	78498.570	0.0000	-0.0665	0.9470	3.602	3.602
29	CREATININE	10	207820.93	0.0000	-203.0302	0.0000	0.536	1.747
30	CREATININE_CMOVAV	10	233789.80	0.0000	-231.1729	0.0000	0.530	1.821
31	CREATININE_CMOVSTD	10	263725.30	0.0000	-213.6878	0.0000	0.110	0.545
32	CVP	10	121734.85	0.0000	-99.7933	0.0000	18.407	44.714
33	CVP_CMOVAV	10	137075.09	0.0000	-116.0867	0.0000	18.119	47.277
34	CVP_CMOVSTD	10	47521.818	0.0000	-73.8026	0.0000	6.298	13.571
35	DERIV1_ABSOLUTE_NEUT_CT	10	63174.243	0.0000	12.0163	0.0000	0.000	0.000
36	DERIV1_ABSOLUTE_NEUT_CT_CMOVAV	10	48955.497	0.0000	58.8763	0.0000	0.000	0.000
37	DERIV1_ABSOLUTE_NEUT_CT_CMOVSTD	10	92626.641	0.0000	88.6559	0.0000	0.000	0.000
38	DERIV1_ALT	10	76121.510	0.0000	37.4061	0.0000	0.000	-0.000
39	DERIV1_ALT_CMOVAV	10	116295.81	0.0000	-125.2304	0.0000	0.000	0.013
40	DERIV1_ALT_CMOVSTD	10	141896.62	0.0000	-121.4974	0.0000	0.003	0.054
41	DERIV1_APTT	9	123445.00	0.0000	13.7649	0.0000	0.000	-0.000

FIG. 10A

12/32

QTMS 4.0 SINGLE FACTOR SCAN SUMMARY TABLE OF PROJECT TESTRUN THE VARIABLES ARE LISTED IN ALPHABETICAL ORDER								
#	VARIABLE	NO. OF INTERVALS	CHI- SQUARE	CHI-SQ. PROB.	T. TEST	T. TEST PROB.	MEAN OF NON.RESP.	MEAN OF RESPONDERS
452	PLATELET_COUNT_CMOVAV	10	202640.10	0.0000	294.3285	0.0000	342.245	169.428
453	PLATELET_COUNT_CMOVSTD	10	119631.52	0.0000	256.7107	0.0000	114.839	60.510
454	PNN1000	3	5189.494	0.0000	65.8434	0.0000	0.014	0.006
455	PNN1000P	4	6197.004	0.0000	25.1068	0.0000	0.024	0.017
456	PNN200	2	20604.205	0.0000	-40.7510	0.0000	0.001	0.007
457	PNN300	2	30763.744	0.0000	-114.6616	0.0000	0.005	0.108
458	PNN400	5	21781.124	0.0000	-120.1509	0.0000	0.139	0.303
459	PNN500	7	8773.508	0.0000	11.5811	0.0000	0.357	0.341
460	PNN600	7	10369.983	0.0000	62.5211	0.0000	0.227	0.162
461	PNN700	6	21305.207	0.0000	215.9426	0.0000	0.122	0.023
462	PNN800	5	15442.091	0.0000	155.1976	0.0000	0.063	0.015
463	PNN900	4	12129.250	0.0000	133.1152	0.0000	0.033	0.008
464	PROTEIN_TOTAL	8	72238.869	0.0000	-108.0584	0.0000	5.147	5.832
465	PROTEIN_TOTAL_CMOVAV	8	67751.460	0.0000	-45.6411	0.0000	5.465	5.678
466	PROTEIN_TOTAL_CMOVSTD	8	169068.31	0.0000	304.3326	0.0000	0.827	0.439
467	RATIO	10	81516.074	0.0000	73.6169	0.0000	0.865	0.548
468	RATIO_CMOVAV	10	144548.40	0.0000	161.7828	0.0000	0.865	0.548
469	RATIO_CMOVSTD	10	46899.538	0.0000	45.6680	0.0000	0.385	0.212
470	RESP_RATE	10	5679.399	0.0000	-35.8331	0.0000	29.038	31.892
471	RESP_RATE_CMOVAV	10	5964.430	0.0000	-40.9464	0.0000	28.881	31.702
472	RESP_RATE_CMOVSTD	10	5824.999	0.0000	-11.7409	0.0000	8.015	8.414
473	SP02	9	84339.963	0.0000	234.2704	0.0000	96.474	89.999
474	SP02_CMOVAV	10	103215.50	0.0000	261.4152	0.0000	96.594	90.052
475	SP02_CMOVSTD	10	5703.834	0.0000	-17.9558	0.0000	2.272	2.436
476	SP02_PULSE_RATE	10	35957.376	0.0000	-208.0506	0.0000	119.222	136.360
477	SP02_PULSE_RATE_CMOVAV	10	44546.210	0.0000	-226.3475	0.0000	119.356	136.539
478	SP02_PULSE_RATE_CMOVSTD	10	64473.674	0.0000	175.6956	0.0000	10.344	5.693
479	TEMPERATURE	9	31893.416	0.0000	62.0223	0.0000	37.184	36.052
480	TEMPERATURE_CMOVAV	10	34350.914	0.0000	58.2470	0.0000	37.224	36.682
481	TEMPERATURE_CMOVSTD	10	26665.736	0.0000	-59.3761	0.0000	0.298	1.000
482	TROPONIN_T	5	1190.672	0.0000	-673.6947	0.0000	0.328	3.955
483	TROPONIN_T_CMOVAV	5	1190.672	0.0000	-1638.000	0.0000	0.303	4.437
484	TROPONIN_T_CMOVSTD	4	1183.480	0.0000	-175.9527	0.0000	0.099	0.683
485	ULF	10	96956.433	0.0000	6.6616	0.0000	0.000	0.000
486	UREA_NITROGEN	9	224418.18	0.0000	-245.6581	0.0000	13.822	49.128
487	UREA_NITROGEN_CMOVAV	10	230750.41	0.0000	-265.8944	0.0000	13.755	50.140
488	UREA_NITROGEN_CMOVSTD	10	250838.84	0.0000	-300.7299	0.0000	5.925	14.020
489	URIC_ACID	5	61063.624	0.0000	497.8559	0.0000	3.985	1.097
490	URIC_ACID_CMOVAV	5	48697.756	0.0000	220.6708	0.0000	3.519	2.798
491	URIC_ACID_CMOVSTD	5	54898.675	0.0000	-132.1616	0.0000	1.063	1.486
492	VHF	10	20886.701	0.0000	-9.9617	0.0000	0.005	0.005

FIG. 10B

13/32

QTMS 4.0	INTERVALS AND THEIR CODED VARIABLES OF PROJECT TESTRUN THAT HAVE -MORTAL- RATES SIGNIFICANTLY HIGHER (WITH 95% CONFIDENCE) THAN THE OVERAL -MORTAL- RATE + 20%						
NEW, CODED VARIABLE	ORIGINAL VARIABLE	INTERVAL	NO. OF PATIENTS-	NO. OF DEADS.	-MORAL- RATE	CONF. INDEX	CONF. INDEX
CODED1	UREA_NITROGEN_CMOVSTD	16.39657339 - 27.47294123	35757	35757	100.00	721	721
CODED2	GGT_CMOVAV	243.375 - 386	19490	19490	100.00	721	721
CODED3	DERIV1_FIBRINOGEN_F_I_CMOVAV	-0.00001388 - 0.000138223	15968	15968	100.00	721	721
CODED4	D_DIMER_QUANT_CMOVSTD	2070.603068 - 6380.224487	15234	15234	100.00	721	721
CODED5	BILIRUBIN_TOTAL_CMOVSTD	1.935597052 - 4.652239604	13845	13845	100.00	721	721
CODED6	DERIV1_FIBRINOGEN_F_I_CMOVSTD	0	13565	13565	100.00	721	721
CODED7	DERIV2_FIBRINOGEN_F_I_CMOVAV	0	13565	13565	100.00	721	721
CODED8	D_DIMER_QUANT_CMOVAV	5323.5	13565	13565	100.00	721	721
CODED9	DERIV2_GGT_CMOVSTD	7.806672E-9 - 1.147579E-8	13133	13133	100.00	721	721
CODED10	DERIV1_ALT_CMOVSTD	0.169115568 - 0.188612207	11552	11552	100.00	721	721
CODED11	DERIV1_BILIRUBIN_TOTAL_CMOVSTD	0.000476013 - 0.000508497	9479	9479	100.00	721	721
CODED12	DERIV1_ALT_CMOVAV	0.039836229 - 0.047007196	9007	9007	100.00	721	721
CODED13	DERIV2_D_DIMER_QUANT_CMOVSTD	8.892813E-7 - 1.061478E-6	7051	7051	100.00	721	721
CODED14	DERIV1_AST_CMOVAV	0.0375103 - 0.041215289	7017	7017	100.00	721	721
CODED15	DERIV1_AST_CMOVSTD	0.15921377 - 0.166908355	7017	7017	100.00	721	721
CODED16	DERIV1_URIC_ACID	0.343344E-7	7017	7017	100.00	721	721
CODED17	DERIV1_URIC_ACID_CMOVSTD	3.990617E-6	7017	7017	100.00	721	721
CODED18	DERIV2_URIC_ACID_CMOVAV	4.79692E-12	7017	7017	100.00	721	721
CODED19	DERIV2_URIC_ACID_CMOVSTD	1.01064E-11	7017	7017	100.00	721	721
CODED20	URIC_ACID	0.9	7017	7017	100.00	721	721
CODED21	DERIV1_FIBRINOGEN_F_I_CMOVAV	-0.00010762 - -0.0000397	6037	6037	100.00	721	721
CODED22	DERIV1_FIBRINOGEN_F_I_CMOVSTD	0.000764479 - 0.000781514	5979	5979	100.00	721	721
CODED23	D_DIMER_QUANT_CMOVSTD	1026.01194 - 2004.861061	5293	5293	100.00	721	721
CODED24	AST_CMOVAV	589.1818182 - 3027	4392	4392	100.00	721	721
CODED25	DERIV1_D_DIMER_QUANT_CMOVAV	0.006040644	4160	4160	100.00	721	721
CODED26	DERIV1_D_DIMER_QUANT_CMOVSTD	0.035055847	4160	4160	100.00	721	721
CODED27	DERIV2_FIBRINOGEN_F_I_CMOVSTD	3.694288E-8	4160	4160	100.00	721	721
CODED28	FIBRINOGEN_F_I_CMOVAV	176.7142857	4160	4160	100.00	721	721
CODED29	DERIV2_FIBRINOGEN_F_I	2.71389E-10 - 4.955707E-9	3963	3963	100.00	721	721
CODED30	DERIV2_D_DIMER_QUANT_CMOVAV	-1.98988E-7 - 4.34546E-8	3697	3697	100.00	721	721
CODED31	D_DIMER_QUANT_CMOVAV	4069.692308 - 4967	3402	3402	100.00	721	721
CODED32	DERIV2_FIBRINOGEN_F_I_CMOVAV	1.616511E-9 - 4.955707E-9	3209	3209	100.00	721	721
CODED33	DERIV1_FIBRINOGEN_F_I	0.000150946 - 0.00043882	3148	3148	100.00	721	721
CODED34	ABSOLUTE_NEUT_CT_CMOVAV	22.125 - 22.6	2910	2910	100.00	721	721
CODED35	DERIV1_BILIRUBIN_TOTAL_CMOVAV	0.000119513 - 0.000125036	2760	2760	100.00	721	721
CODED36	DERIV2_BILIRUBIN_TOTAL_CMOVSTD	9.679595E-9 - 9.902011E-9	2760	2760	100.00	721	721
CODED37	AST	2811 - 2860	2195	2195	100.00	721	721
CODED38	DERIV2_AST_CMOVSTD	3.286978E-6	1265	1265	100.00	721	721
CODED39	DERIV1_TROPONIN_T	0.000039239	207	207	100.00	721	721
CODED40	DERIV1_TROPONIN_T_CMOVAV	0.000039239	207	207	100.00	721	721
CODED41	DERIV1_TROPONIN_T_CMOVSTD	0	207	207	100.00	721	721
CODED42	DERIV2_TROPONIN_T	0	207	207	100.00	721	721
CODED43	DERIV2_TROPONIN_T_CMOVAV	0	207	207	100.00	721	721

FIG. 11A

14/32


QTMS 4.0	INTERVALS AND THEIR CODED VARIABLES OF PROJECT TESTRUN THAT HAVE -MORTAL- RATES SIGNIFICANTLY HIGHER (WITH 95% CONFIDENCE) THAN THE OVERAL -MORTAL- RATE + 20%						
NEW, CODED VARIABLE	ORIGINAL VARIABLE	INTERVAL	NO. OF PATIENTS-	NO. OF DEADS.	-MORTAL- RATE	CONF. INDEX	CONF. INDEX
CODED2366	DERIV2_MULTI_HEMA	-8.0141E-10 - -2.9432E-10	39920	6915	17.32	125	122
CODED2367	RESP_RATE	47.6779661 - 199.1208054	69841	12026	17.22	124	122
CODED2368	DERIV2_MULTI_PC02_VEN_CMOVAV	-7.60128E-9 - -2.9522E-9	23376	4073	17.42	126	122
CODED2369	DERIV2_UREA_NITROGEN_CMOVSTD	5.88583E-10 - 9.20105E-10	40019	6925	17.30	125	122
CODED2370	MULTI_BASE_CMOVSTD	4.122256279 - 4.555373593	33774	5853	17.33	125	122
CODED2371	DERIV2_PHOSPHORUS	-1.5461E-10 - -6.4608E-11	39303	6796	17.29	125	122
CODED2372	DERIV2_PLATELET_COUNT_CMOVSTD	9.06034E-8 - 1.924971E-7	38902	6722	17.28	125	122
CODED2373	TEMPERATURE	38.50048732 - 73.48275862	10510	1853	17.63	127	122
CODED2374	DERIV1_MULTI_PO2_VEN	-0.00006551 - -1.12177E-6	23935	4156	17.36	125	122
CODED2375	MULTI_PC02_VEN_CMOVSTD	10.95445115 - 13.83032096	23482	4078	17.37	125	122
CODED2376	PHOSPHORUS	5.7 - 6.1	46887	8071	17.21	124	122
CODED2377	MULTI_PC02	18 - 37	25502	4421	17.34	125	122
CODED2378	DERIV1_MULTI_ART_CMOVAV	0.010179604 - 0.490510929	23861	4140	17.35	125	122
CODED2379	NIBP_SYS_CMOVAV	81.00053763 - 88.13232323	71684	12285	17.14	124	122
CODED2380	DERIV1_RATIO_CMOVAV	0.0000E-309 - 9.001911E-6	76232	13058	17.13	124	122
CODED2381	SP02_PULSE_RATE_CMOVAV	141.7871345 - 151.7447355	75290	12882	17.11	123	121
CODED2382	DERIV1_MULTI_PH_CMOVAV	5.026518E-7 - 1.763157E-6	25885	4478	17.30	125	121
CODED2383	DERIV1_CVP_CMOVSTD	0.02994997 - 0.040792	34016	5860	17.23	124	121
CODED2384	MULTI_PC02_CMOVSTD	6.488084316 - 7.789054225	31857	5491	17.24	124	121
CODED2385	MULTI_I0N_CMOVSTD	0.714169555 - 1.814873399	25594	4412	17.24	124	121
CODED2386	DERIV2_MULTI_POT_CMOVAV	9.227356E9 - 1.010527E-7	39964	6851	17.14	124	121
CODED2387	DERIV2_MULTI_BICARB_CMOVSTD	1.692021E-9 - 6.077E-7	33410	5736	17.17	124	121
CODED2388	DERIV2_MULTI_BASE	4.69537E10 - 1.913356E-9	33769	5793	17.15	124	121
CODED2389	NN_INTERVAL	0.391617667 - 0.424661198	84693	14385	16.98	123	121
CODED2390	MULTI_PO2_VEN_CMOVAV	44.16 - 45.475	23473	4039	17.21	124	121
CODED2391	DERIV2_CHLORIDE	-1.46469E-9 - -7.4057E-10	39332	6715	17.07	123	120
CODED2392	DERIV1_MULTI_I0N	-4.59294E-6 - -7.7983E-7	32487	5558	17.11	123	120
CODED2393	DERIV2_RATIO	-0.00001545 - -6.32576E-6	76232	12928	16.96	122	120
CODED2394	MULTI_PO2	115 - 146	24239	4161	17.17	124	120
CODED2395	O2_SATURATON_CMOVSTD	3.770735608 - 4.706207905	28442	4872	17.13	124	120
CODED2396	DERIV1_CREATININE	-7.76159E-7 - -3.64093E-7	39410	6718	17.05	123	120
CODED2397	HIBP_MEAN_CMOVAV	69.00014798 - 74	72791	12334	16.94	122	120
CODED2398	DERIV2_MULTI_POT_CMOVSTD	2.047386E-8 - 3.915641E-8	40006	6808	17.02	123	120
CODED2399	CALCIUM	9.4 - 9.6	39741	6759	17.01	123	120
CODED2400	PHOSPHORUS_CMOVSTD	0.98386991 - 1.063810176	38597	6566	17.01	123	120

FIG. 11B



MULTI_SODIUM RAW OPTIMIZED-EMITTING

Int	RANGE	PICU COUNT	DEAD COUNT	MORTALITY PERCENT	ALARM
1	101 - 132	734	240	30.61	FLAT
2	133 - 134	802	162	20.20	LOW
3	135 - 136	1,338	192	14.35	LOW
4	131 - 138	1,365	190	13.92	LOW
5	139 - 140	1,042	194	10.02	LOW
6	141 - 143	935	230	26.60	FLAT
7	144 - 149	756	366	47.78	HIGH
8	150 - 177	530	372	76.40	HIGH

FIG. 12A



15/32



16/32

POWER PLAY WEB EXPLORER

DEADDAY

PATIENT

ORGAN

RENAL

RAW

MOVSTD

UREA_NITROGEN_CMOVSTD

"UREA_NITROGEN_CMOVSTD: 16.39657339 - 27.47294123"

"UREA_NITROGEN_CMOVSTD: 3.342709617 - 3.933264586"

"UREA_NITROGEN_CMOVSTD: 8.534523726 - 10.59918879"

"UREA_NITROGEN_CMOVSTD: 10.63427015 - 16.32218877"

MULTI_ION_CMOVSTD

"MULTI_ION_CMOVSTD: 0.418296519 - 0.489322257"

"MULTI_ION_CMOVSTD: 0.225907523 - 0.248003889"

"MULTI_ION_CMOVSTD: 0.292217229 - 0.350749054"

"MULTI_ION_CMOVSTD: 0.714169555 - 1.614673399"

"MULTI_ION_CMOVSTD: 0.35097879 - 0.418294155"

"MULTI_ION_CMOVSTD: 0.489621409 - 0.714144674"

URIC_ACID_CMOVSTD

CHLORIDE_CMOVSTD

CALCIUM_CMOVSTD

MULTI_SODIUM_CMOVSTD

PHOSPHORUS_CMOVSTD

MULTI_GLUCOSE_CMOVSTD

CREATININE_CMOVSTD

MULTI_POT_CMOVSTD

RAW

MULTI_POT

MAGNESIUM

PHOSPHORUS

CALCIUM

CREATININE

MULTI_ION

CO2

MULTI_GLUCOSE

URIC_ACID

MULTI_SODIUM

CHLORIDE

UREA_NITROGEN

MOVAVE

MULTI_SODIUM_CMOVAV

CHLORIDE_CMOVAV

CALCIUM_CMOVAV

UREA_NITROGEN_CMOVAV

MULTI_GLUCOSE_CMOVAV

URIC_ACID_CMOVAV

MULTI_POT_CMOVAV

MAGNESIUM_CMOVAV

CO2_CMOVAV

MULTI_ION_CMOVAV

PHOSPHORUS_CMOVAV

CREATININE_CMOVAV

DERV1

DERV2

DEADDAY

PATIENT

ORGAN

PATIENT DATE

VAR_TYPE

VAR_SUBTYPE

VARIABLE

LABEL

RENAL

HEPATIC

CARDIAC

RESPIRATOR

HEMATOLOGY

ID

uv 17a8d942

RAW

RAW

182,890

29,670

425,297

193,163

87,295

33,778

MOVAVE

182,905

27,946

141,507

180,736

94,055

35,706

MOVSTD

155,469

29,686

92,231

174,997

93,122

34,270

VAR_SUBTYPE

521,264

87,302

659,035

548,896

274,472

103,754

DERV1

RAW

157,004

29,127

122,000

154,691

105,452

53,148

MOVAVE

158,898

29,698

99,663

154,701

108,209

31,757

MOVSTD

162,788

29,062

118,837

170,037

104,810

18,485

VAR_SUBTYPE

478,690

87,887

340,600

499,429

318,471

103,390

DERV2

RAW

170,485

38,553

90,961

176,898

91,401

49,880

MOVAVE

168,423

41,675

97,201

171,856

92,949

32,054

MOVSTD

174,876

33,040

124,661

181,030

88,940

18,466

VAR_SUBTYPE

513,784

113,268

312,823

529,764

273,230

100,400

VAR_TYPE

1,513,738

288,457

1,312,358

1,578,109

866,233

307,544

uv 19732a73

RAW

RAW

47,301

0

64,635

67,941

11,705

11,231

MOVAVE

45,553

0

21,386

59,110

13,310

12,711

MOVSTD

41,383

0

14,979

59,378

11,302

9,643

VAR_SUBTYPE

134,237

0

104,000

186,429

36,317

33,585

DERV1

RAW

51,438

0

26,688

59,503

26,036

13,057

MOVAVE

57,243

0

22,000

53,249

23,553

12,988

MOVSTD

54,061

0

18,708

59,812

26,275

12,957

VAR_SUBTYPE

162,742

0

67,404

172,564

75,864

39,002

DERV2

RAW

60,777

0

19,181

64,667

19,197

14,778

MOVAVE

53,898

0

20,892

65,511

19,571

12,957

MOVSTD

47,654

0

19,385

57,749

19,571

12,807

VAR_SUBTYPE

162,329

0

59,761

187,927

59,339

40,542

VAR_TYPE

459,308

0

231,165

546,920

170,520

113,129

uv 3b9967bb

RAW

RAW

41,418

6,722

112,551

42,203

12,080

3,335

MOVAVE

40,233

7,780

33,361

41,402

13,376

3,611

MOVSTD

32,662

9,031

37,226

42,190

9,572

2,430

VAR_SUBTYPE

114,313

23,533

183,138

125,795

35,028

9,376

DERV1

RAW

38,507

6,702

42,925

44,291

17,940

5,745

MOVAVE

42,565

7,428

36,929

43,953

17,179

4,762

MOVSTD

30,632

5,835

43,529

44,276

11,783

724

VAR_SUBTYPE

111,704

19,965

129,283

132,523

46,902

11,231

DERV2

RAW

42,454

9,031

37,086

46,383

11,996

4,763

MOVAVE

36,076

5,835

32,595

48,434

7,952

4,763

MOVSTD

38,924

8,337

46,344

49,991

7,639

2,163

VAR_SUBTYPE

117,454

23,203

116,025

144,818

27,587

11,689

VAR_TYPE

313,471

66,701

422,446

403,131

109,517

32,296

uv 3f700bet

RAW

RAW

135,073

1,834

159,643

123,972

96,196

27,568

MOVAVE

149,924

1,834

77,910

121,157

97,078

31,329

MOVSTD

107,264

1,034

71,000

107,720

84,112

20,100

FIG. 14A

17/32

POWER PLAY WEB EXPLORER

NEWLIVE

PATIENT

ORGAN

RENAL

RAW

MOVSTD

UREA_NITROGEN_CMOVSTD

"UREA_NITROGEN_CMOVSTD: 16.39657339 - 27.47294123"

"UREA_NITROGEN_CMOVSTD: 3.342789617 - 3.983264586"

"UREA_NITROGEN_CMOVSTD: 8.534523726 - 10.59918879"

"UREA_NITROGEN_CMOVSTD: 10.63427015 - 16.32218877"

MULTI_ION_CMOVSTD

"MULTI_ION_CMOVSTD: 0.418296519 - 0.489322257"

"MULTI_ION_CMOVSTD: 0.225907523 - 0.248003889"

"MULTI_ION_CMOVSTD: 0.292217229 - 0.350749054"

"MULTI_ION_CMOVSTD: 0.714169555 - 1.614673399"

"MULTI_ION_CMOVSTD: 0.35097879 - 0.418294155"

"MULTI_ION_CMOVSTD: 0.489621409 - 0.714144674"

URIC_ACID_CMOVSTD

"URIC ACID_CMOVSTD: 0.892728562 - 1.4591424"

"URIC ACID_CMOVSTD: 1.5264398937"

CHLORIDE_CMOVSTD

"CHLORIDE_CMOVSTD: 2.515148907 - 3.082207001"

"CHLORIDE_CMOVSTD: 4.140393356 - 4.581805087"

"CHLORIDE_CMOVSTD: 5.642149072 - 6.358346596"

"CHLORIDE_CMOVSTD: 4.502575695 - 4.97214463"

"CHLORIDE_CMOVSTD: 6.882104574 - 16.10279479"

"CHLORIDE_CMOVSTD: 6.361341578 - 6.858296573"

"CHLORIDE_CMOVSTD: 4.97352297 - 5.631988836"

CALCIUM_CMOVSTD

"CALCIUM_CMOVSTD: 0.83597676 - 1.117735807"

"CALCIUM_CMOVSTD: 0 - 0.357628974"

"CALCIUM_CMOVSTD: 0.262769136 - 0.4"

"CALCIUM_CMOVSTD: 0.697186374 - 0.83522676"

"CALCIUM_CMOVSTD: 0.588816924 - 0.645638604"

MULTI_SODIUM_CMOVSTD

"MULTI_SODIUM_CMOVSTD: 5.100238169 - 5.884365159"

"MULTI_SODIUM_CMOVSTD: 0 - 2.004540301"

"MULTI_SODIUM_CMOVSTD: 3.332768311 - 3.879432948"

"MULTI_SODIUM_CMOVSTD: 3.880415844 - 4.343167405"

"MULTI_SODIUM_CMOVSTD: 5.890259476 - 12.81670227"

"MULTI_SODIUM_CMOVSTD: 4.345321413 - 5.09835467"

PHOSPHORUS_CMOVSTD

"PHOSPHORUS_CMOVSTD: 1.423345355 - 1.84327075"

"PHOSPHORUS_CMOVSTD: 0.843528137 - 0.982542422"

"PHOSPHORUS_CMOVSTD: 1.064489614 - 1.219289411"

"PHOSPHORUS_CMOVSTD: 0.98386991 - 1.063810176"

"PHOSPHORUS_CMOVSTD: 1.85265789 - 2.969848481"

"PHOSPHORUS_CMOVSTD: 1.220584115 - 1.422176014"

"PHOSPHORUS_CMOVSTD: 0.7288058 - 0.847068411"

MULTI_GLUCCOSECMOVSTD

CREATININE_CMOVSTD

MULTI_POT_CMOVSTD

RAW

PATIENT - ORGAN - PATIENT DATE - VAR_TYPE - VAR_SUBTYPE - VARIABLE - LABEL

EVENT AS VALUES

RENAL

HEPATIC

CARDIAC

RESPIRATOR

HEMATOLOGY

ID

ORGAN

uv 0c43e928

RAW

RAW

0

0

89,955

2,971

0

0

92,926

MOVAVE

0

0

27,633

2,849

0

0

30,482

MOVSTD

0

0

11,017

3,996

0

0

15,013

VAR_SUBTYPE

0

0

128,605

9,816

0

0

138,421

DERV1

RAW

0

0

24,907

1,792

0

0

26,699

MOVAVE

0

0

17,424

784

0

0

18,208

MOVSTD

0

0

17,361

1,501

0

0

18,662

VAR_SUBTYPE

0

0

59,692

4,077

0

0

63,769

DERV2

RAW

0

0

14,123

1,749

0

0

15,872

MOVAVE

0

0

17,260

704

0

0

18,044

MOVSTD

0

0

19,245

1,500

0

0

20,745

VAR_SUBTYPE

0

0

50,628

4,033

0

0

54,661

VAR_TYPE

0

0

238,925

17,926

0

0

256,851

uv 19cd239b

RAW

RAW

69,114

0

53,916

20,413

5,956

5,350

154,749

MOVAVE

50,948

0

4,681

1,773

7,251

9,428

74,081

MOVSTD

35,848

0

4,298

19,206

10,226

4,731

74,309

VAR_SUBTYPE

155,910

0

62,095

41,392

23,483

19,509

303,139

DERV1

RAW

75,262

0

25,760

23,834

18,106

8,741

151,703

MOVAVE

40,131

0

12,029

34,343

13,722

7,623

107,848

MOVSTD

29,768

0

16,578

11,162

9,403

2,135

68,066

VAR_SUBTYPE

145,181

0

53,367

69,339

41,231

18,499

327,617

DERV2

RAW

73,889

0

12,195

26,802

20,298

8,659

141,843

MOVAVE

79,191

0

15,140

19,521

4,390

3,261

121,503

MOVSTD

43,518

0

16,768

14,929

1,805

9,056

86,076

VAR_SUBTYPE

196,598

0

44,103

61,252

26,493

20,976

349,422

VAR_TYPE

497,698

0

150,365

171,983

91,157

58,984

980,178

uv 352ed5e

RAW

RAW

57,320

0

97,582

67,287

36,580

0

258,769

MOVAVE

67,254

0

29,340

71,282

38,212

0

205,108

MOVSTD

58,307

0

15,308

52,930

40,421

4,483

171,449

VAR_SUBTYPE

182,881

0

141,230

191,499

115,233

4,489

635,326

DERV1

RAW

75,158

0

30,369

55,738

51,485

14,633

227,383

MOVAVE

42,465

0

36,860

47,792

27,463

9,558

164,138

MOVSTD

68,263

0

21,303

38,234

41,138

4,483

173,421

VAR_SUBTYPE

185,886

0

88,532

141,764

120,086

28,674

664,942

DERV2

RAW

66,403

0

13,742

49,325

38,748

8,843

177,061

MOVAVE

46,998

0

22,143

60,902

33,120

10,380

173,543

MOVSTD

55,409

0

22,688

47,792

20,473

9,558

155,920

VAR_SUBTYPE

168,810

0

58,573

158,019

92,341

28,781

506,524

VAR_TYPE

537,577

0

288,335

491,282

327,660

61,938

1,706,792

uv 28ba5e84

RAW

RAW

67,871

10,622

54,549

63,027

14,420

15,171

225,660

MOVAVE

38,904

10,649

35,594

31,889

0

38,888

FIG. 14B


18/32

PATIENT	STARTDATE	P_DEAD1	ID	CARDIAC	RENAL	RESPIRATORY,	HEPATIE	HEMATOLOGY
uv_17a8d942	17JUL	0.837234	0.0	40.0	0.0	18.0	0.0	0.0
uv_17a8d942	18JUL	0.818931	0.0	40.2	0.0	18.0	0.0	0.0
uv_17a8d942	19JUL	0.831449	0.0	37.5	0.0	18.0	0.0	0.0
uv_17a8d942	20JUL	0.660144	0.0	37.8	0.0	18.0	0.0	0.0
uv_17a8d942	21JUL	0.720858	0.0	38.9	0.0	18.0	0.0	0.0
uv_17a8d942	22JUL	0.645172	0.0	37.1	0.0	18.0	0.0	0.0
uv_17a8d942	23JUL	0.746791	0.0	36.4	0.0	18.0	0.0	0.0
uv_17a8d942	24JUL	0.869837	0.0	41.4	5.3	8.4	0.0	0.0
uv_17a8d942	25JUL	0.998454	7.6	43.7	67.2	30.4	0.0	32.4
uv_17a8d942	26JUL	0.999996	18.5	48.6	79.1	39.6	0.0	42.9
uv_17a8d942	27JUL	0.999999	17.2	49.2	82.2	37.6	8.1	42.4
uv_17a8d942	28JUL	1.000000	16.2	52.0	84.7	39.2	10.0	43.7
uv_17a8d942	29JUL	1.000000	16.0	51.4	84.2	40.6	10.0	44.8
uv_17a8d942	30JUL	1.000000	16.8	53.0	83.2	41.4	10.0	46.6
uv_17a8d942	31JUL	1.000000	16.2	51.3	86.8	39.7	10.0	44.1
uv_17a8d942	01AUG	1.000000	18.8	55.9	85.5	54.5	10.0	46.4
uv_17a8d942	02AUG	1.000000	17.7	56.4	85.1	70.3	10.0	46.5
uv_17a8d942	03AUG	1.000000	17.1	62.4	82.7	103.3	10.0	46.0
uv_17a8d942	04AUG	1.000000	17.6	60.9	85.3	105.7	10.0	46.6
uv_17a8d942	05AUG	0.999961	18.0	62.9	85.8	110.3	10.0	45.4
uv_17a8d942	06AUG	0.999848	17.2	66.4	85.1	105.0	10.0	46.6
uv_17a8d942	07AUG	0.999999	16.1	71.6	86.5	103.4	16.6	46.3
uv_17a8d942	08AUG	1.000000	17.0	68.8	86.0	103.9	35.9	45.2
uv_17a8d942	09AUG	1.000000	16.7	69.1	85.5	98.4	43.8	59.5
uv_17a8d942	10AUG	1.000000	17.8	70.3	83.9	99.7	44.9	65.7
uv_17a8d942	11AUG	1.000000	16.0	67.4	87.4	101.3	44.6	65.1
uv_17a8d942	12AUG	1.000000	18.0	64.7	84.7	105.4	36.1	65.5
uv_17a8d942	13AUG	1.000000	17.7	66.5	83.8	103.5	23.3	60.9
uv_17a8d942	14AUG	1.000000	17.4	65.8	83.4	101.1	12.9	44.7
uv_17a8d942	15AUG	1.000000	17.3	62.6	81.5	101.8	10.0	43.7
uv_17a8d942	16AUG	1.000000	17.0	65.0	81.1	103.3	10.0	42.5
uv_17a8d942	17AUG	1.000000	17.0	64.3	77.7	102.8	10.0	41.1
uv_17a8d942	18AUG	1.000000	17.3	60.0	77.6	86.9	4.5	36.7
uv_17a8d942	19AUG	0.925828	19.0	52.2	40.1	56.7	0.0	23.9

FIG. 14C

19/32

Int	RANGE	PICU COUNT	DEAD COUNT	MORTALITY PERCENT	ALARM
1	123.25 - 154.7372737	36,386	7,707	21.18	FLAT
2	124.7333323 - 135.7666037	36,438	10,775	29.55	HIGH
3	133.7777378 - 136.0452331	36,400	10,477	35.70	HIGH
4	130.0967742 - 136.04666657	36,389	6,384	17.47	LOW
5	138.6727253 - 137.0633333	28,592	977	2.90	LOW
6	137.0009058 - 137.6638462	36,360	1,795	24.96	LOW
7	1370004667 - 135	37,547	2,719	24.96	LOW
8	1390046667 - 100.3444478	30,404	2,877	7.99	LOW
9	140.30 - 144.0287190	75,950	7,470	30.42	FLAT
10	144083377 - 159	30,417	38,334	30.47	HIGH

FIG. 15A**FIG. 15B**

20/32

Int	RANGE	PICU COUNT	DEAD COUNT	MORTALITY PERCENT	ALARM
1	0 1.000207398	36,150	11,500	12.18	HIGH
2	1.20000800 - 2300401077	36,204	210	0.27	LOW
3	331.3521322 - 2030370312	36,148	408	1.03	LOW
4	2344601122 - 3099707539	36,003	1,419	3.24	LOW
5	3.10166713 - 3.642343568	35,909	7,138	19.81	FLAT
6	3.545715605 - 4083188015	36,3467	5,254	14.57	LOW
7	4.4859826100 - 4.490475027	36,007	13,006	30.26	HIGH
8	4.498720102 - 5.25237312	36,003	10,302	29.17	HIGH
9	4.498720103 - 5.938303117	36,383	14,575	40.00	HIGH
10	5.939778848 - 12.88670227	44,926	19,002	56.28	HIGH

FIG. 15C




FIG. 15D




FIG. 16A

22/32




FIG. 16B

23/32




FIG. 17A

24/32





FIG. 17B

25/32



<?xml version="1.0" encoding="utf-8" standalone="yes"?>

<PMML version="3.1">

<xmlns="http://www.dmg.org/PMML-3.1">

<xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance">

<Header copyright="copyright(c) 2002 SAS INSTITUTE Inc., Cary, NC, USA> All Rights Reserved.">

<Application name="SAS(r)" version="9.1"/>

<Timestamp>2009-07-09 10:28:15</Timestamp>

</Header>

<DataDictionary numberOfFields="2401">

<DataField name="CODED1" displayName="UREA NITROGEN CMOVSTD: 16.39657339 - 27.47294123" optype="categorical" datatype="double"/>
 <DataField name="CODED10" displayName="DERIV1 ALT CMOVSTD: 0.169115568 - 0.188912207" optype="categorical" datatype="double"/>
 <DataField name="CODED100" displayName="DERIVT HR ECG 2: -0.00017142 - 6.765583E-7" optype="categorical" datatype="double"/>
 <DataField name="CODED1000" displayName="DERIV2 MULTI PCO2 CMOVSTD: 4.472721E-9 - 5.062997E-8" optype="categorical" datatype="double"/>
 <DataField name="CODED1001" displayName="DERIV2 INR CMOVAV: -4.6668E-11 - 9.5463E-13" optype="categorical" datatype="double"/>
 <DataField name="CODED1002" displayName="DERIV1 CREATININE CMOVAV: 3.042728E-7 - 9.530929E-7" optype="categorical" datatype="double"/>
 <DataField name="CODED1003" displayName="DERIV1 CO2: -0.00005617 - 0.00002828" optype="categorical" datatype="double"/>
 <DataField name="CODED1004" displayName="DERIV1 PROTEIN TOTAL CMOVSTD: 0.00001621 - 0.000229466" optype="categorical" datatype="double"/>
 <DataField name="CODED1005" displayName="DERIV1 MULTI HEMA: -0.02649999 - 0.00011777" optype="categorical" datatype="double"/>
 <DataField name="CODED1006" displayName="PNN400: 0.992217899 - 0.99397504" optype="categorical" datatype="double"/>
 <DataField name="CODED1007" displayName="DERIV1 MULTI BASE: -0.00012791 - 0.00005334" optype="categorical" datatype="double"/>
 <DataField name="CODED1008" displayName="DERIV1 ETCO2: 0.013668688 - 0.033916317" optype="categorical" datatype="double"/>
 <DataField name="CODED1009" displayName="MULTI POT: 4.6 - 5" optype="categorical" datatype="double"/>
 <DataField name="CODED101" displayName="DERIV2 D DIMER QUANT: -8.20058E-7 - 2.66792E-7" optype="categorical" datatype="double"/>
 <DataField name="CODED1010" displayName="DERIV2 ETCO2 CMOVAV: 0.000109508 - 0.256615646" optype="categorical" datatype="double"/>
 <DataField name="CODED1011" displayName="DERIV1 MULTI POT CMOVSTD: 0.000099646 - 0.000126988" optype="categorical" datatype="double"/>
 <DataField name="CODED1012" displayName="DERIV1 MULTI ART CMOVAV: -0.0005293 - 0.000954294" optype="categorical" datatype="double"/>
 <DataField name="CODED1013" displayName="DERIV2 ETCO2: -0.00158687 - 0.00063708" optype="categorical" datatype="double"/>
 <DataField name="CODED1014" displayName="DERIV2 MULTI HEMA CMOVSTD: 1.426777E-8 - 3.857822E-8" optype="categorical" datatype="double"/>
 <DataField name="CODED1015" displayName="DERIV2 ETCO2: 0.000563123 - 0.001238681" optype="categorical" datatype="double"/>
 <DataField name="CODED1016" displayName="DERIV2 MULTI SODIUM: 7.32382E-10 - 6.369898E-9" optype="categorical" datatype="double"/>
 <DataField name="CODED1017" displayName="O2 SATURATION: 99.3 - 99.8" optype="categorical" datatype="double"/>
 <DataField name="CODED1018" displayName="DERIV1 ETCO2 CMOVSTD: 0.021234202 - 0.027710075" optype="categorical" datatype="double"/>
 <DataField name="CODED1019" displayName="DERIV1 MAGNESIUM CMOVSTD: 0.000024436 - 0.000036305" optype="categorical" datatype="double"/>
 <DataField name="CODED102" displayName="DERIV2 HR ECG 2: -0.00001333 - 2.660875E-6" optype="categorical" datatype="double"/>
 <DataField name="CODED1020" displayName="DERIV2 MULTI ART: -0.00077248 - 0.00019346" optype="categorical" datatype="double"/>
 <DataField name="CODED1021" displayName="DERIV2 MULTI TON: 5.40931E-9 - 0.038393959" optype="categorical" datatype="double"/>
 <DataField name="CODED1022" displayName="MULTI HEMA CMOVSTD: 2.621340987 - 3.074556879" optype="categorical" datatype="double"/>
 <DataField name="CODED1023" displayName="DERIV2 ABSOLUTE NEUT CT: -2.31162E-6 - 6.7396E-10" optype="categorical" datatype="double"/>
 <DataField name="CODED1024" displayName="MULTI BASE CMOVSTD: 4.571765489 - 5.274077818" optype="categorical" datatype="double"/>
 <DataField name="CODED1025" displayName="DERIVT ETCO2 CMOVSTD: 0.027710256 - 0.034521434" optype="categorical" datatype="double"/>
 <DataField name="CODED1026" displayName="INR CMOVAV: 0.95 - 1.1873" optype="categorical" datatype="double"/>
 <DataField name="CODED1027" displayName="DERIV1 PLATELET COUNT: 0.000243019 - 0.000487572" optype="categorical" datatype="double"/>
 <DataField name="CODED1028" displayName="DERIV1 MULTI PCO2 CMOVAV: -0.00002564 - 0.000011775" optype="categorical" datatype="double"/>
 <DataField name="CODED1029" displayName="DERIV1 MULTI PO2 CMOVAV: -3.51260854 - 0.00244864" optype="categorical" datatype="double"/>
 <DataField name="CODED103" displayName="FIBRINOGEN F I CMOVSTD: 63.00905511 - 69.99965986" optype="categorical" datatype="double"/>
 <DataField name="CODED1030" displayName="UREA NITROGEN CMOVSTD: 3.342769617 - 3.983264586" optype="categorical" datatype="double"/>
 <DataField name="CODED1031" displayName="DERIV2 CALCIUM: 3.17811E-8 - 2.2886E-10" optype="categorical" datatype="double"/>
 <DataField name="CODED1032" displayName="DERIV1 HEART RATE: -0.00018537 - 0.000171563" optype="categorical" datatype="double"/>
 <DataField name="CODED1033" displayName="MULTI TON CMOVSTD: 0.418296519 - 0.489322257" optype="categorical" datatype="double"/>
 <DataField name="CODED1034" displayName="MULTI HEMO: 12.8 - 13.7" optype="categorical" datatype="double"/>
 <DataField name="CODED1035" displayName="DERIVT PLATELET COUNT CMOVSTD: 0.002272894 - 0.005793821" optype="categorical" datatype="double"/>
 <DataField name="CODED1036" displayName="MAGNESIUM: 2.2 - 2.3" optype="categorical" datatype="double"/>
 <DataField name="CODED1037" displayName="MULTI SODIUM CMOVAV: 135.8571429 - 136.2666667" optype="categorical" datatype="double"/>
 <DataField name="CODED1038" displayName="DERIV2 CO2 CMOVAV: 3.36591E-10 - 9.83044E-10" optype="categorical" datatype="double"/>
 <DataField name="CODED1039" displayName="DERIV2 ETCO2 CMOVSTD: 0.000926837 - 0.00122" optype="categorical" datatype="double"/>
 <DataField name="CODED104" displayName="DERIV2 URIC ACID: 2.97076E-13 - 3.06137E-12" optype="categorical" datatype="double"/>
 <DataField name="CODED1040" displayName="MULTI ART CMOVSTD: 1.662423176 - 2.336502839" optype="categorical" datatype="double"/>
 <DataField name="CODED1041" displayName="DERIVT TEMPERATURE: -0.00129574 - 0.00009641" optype="categorical" datatype="double"/>
 <DataField name="CODED1042" displayName="DERIV2 PROTEIN TOTAL: -2.2461E-12 - 1.5411E-12" optype="categorical" datatype="double"/>
 <DataField name="CODED1043" displayName="DERIV2 PROTEIN TOTAL CMOVAV: 1.41463E-10 - 4.38957E-10" optype="categorical" datatype="double"/>
 <DataField name="CODED1044" displayName="DERIV1 MULTI PH CMOVSTD: 0.000049497 - 0.000241969" optype="categorical" datatype="double"/>
 <DataField name="CODED1045" displayName="DERIV2 MULTI LACTATE CMOVSTD: 2.923545E-6 - 0.000072368" optype="categorical" datatype="double"/>
 <DataField name="CODED1046" displayName="HEART RATE CMOVAV: 153.4790323 - 268.9025424" optype="categorical" datatype="double"/>
 <DataField name="CODED1047" displayName="DERIVT TEMPERATURE CMOVAV: 2.10533E-17 - 0.000028906" optype="categorical" datatype="double"/>
 <DataField name="CODED1048" displayName="TEMPERATURE CMOVAV: 37.55041667 - 37.89931773" optype="categorical" datatype="double"/>
 <DataField name="CODED1049" displayName="PHOSPHORUS: 6.8 - 9.2" optype="categorical" datatype="double"/>
 <DataField name="CODED105" displayName="URIC ACID CMOVSTD: 0.892728562 - 1.45951424" optype="categorical" datatype="double"/>
 <DataField name="CODED1050" displayName="DERIV1 MULTI PH: 2.653449E-6 - 0.000011868" optype="categorical" datatype="double"/>
 <DataField name="CODED1051" displayName="DERIV1 MULTI POT: -0.00007468 - 0.00003677" optype="categorical" datatype="double"/>
 <DataField name="CODED1052" displayName="TEMPERATURE CMOVSTD: -0.00010851 - 0.021978418" optype="categorical" datatype="double"/>

FIG. 18A

27/32

LMPMML


```

<?xml version="1.0" encoding="utf-8" standalone="yes"?>
<PMML version="3.1"
  <xmlns="http://www.dmg.org/PMML-3.1"
  <xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance">
  <Header copyright="copyright(c) 2002 SAS INSTITUTE Inc., Cary, NC, USA> All Rights Reserved.">
  <Application name="SAS(r)" version="9.1"/>
  <Timestamp>2009-07-09 09:56:09</Timestamp>
  </Header>
  <DataDictionary numberFields="1033">
    <DataField name="CODED1" displayName="DERIV2 MULTI ION CMOVAV: 1.734399E-8 - 1.776157E-8" optype="categorical" datatype="double"/>
    <DataField name="CODED10" displayName="DERIVT MULTI LACTATE CMOVAV: -5.56887E-6" optype="categorical" datatype="double"/>
    <DataField name="CODED100" displayName="MULTI LACTATE CMOVSTD: 0.47934547 - 0.480319787" optype="categorical" datatype="double"/>
    <DataField name="CODED1000" displayName="DERIV1 NN INTERVAL CMOVSTD: 0.236848969 - 0.368501314" optype="categorical" datatype="double"/>
    <DataField name="CODED1001" displayName="DERIV1 MULTI LACTATE CMOVAV: 0.000010215 - 0.000099244" optype="categorical" datatype="double"/>
    <DataField name="CODED1002" displayName="MULTI PO2 VEN: 32 - 33" optype="categorical" datatype="double"/>
    <DataField name="CODED1003" displayName="DERIV2 HEART RATE CMOVSTD: 0.006816331 - 0.008036703" optype="categorical" datatype="double"/>
    <DataField name="CODED1004" displayName="HEART RATE CMOVSTD: 3.579963781 - 4.807204951" optype="categorical" datatype="double"/>
    <DataField name="CODED1005" displayName="DERIVT SPO2 PULSE RATE: -0.12506479 - -0.05671283" optype="categorical" datatype="double"/>
    <DataField name="CODED1006" displayName="DERIV2 SPO2 PULSE RATE: 0.004817804 - 0.011495646" optype="categorical" datatype="double"/>
    <DataField name="CODED1007" displayName="DERIV2 HR ECG 2: -0.00879819 - -0.00444886" optype="categorical" datatype="double"/>
    <DataField name="CODED1008" displayName="RATIO: 0.377831205 - 0.456531693" optype="categorical" datatype="double"/>
    <DataField name="CODED1009" displayName="NIBP MEAN: 72.00103684 - 77" optype="categorical" datatype="double"/>
    <DataField name="CODED101" displayName="MULTI PCO2 VEN CMOVAV: 71.2173913 - 72.31818182" optype="categorical" datatype="double"/>
    <DataField name="CODED1010" displayName="DERIV1 MULTI ION CMOVSTD: 0.000407608 - 0.000700346" optype="categorical" datatype="double"/>
    <DataField name="CODED1011" displayName="DERIV2 RESP RATE CMOVSTD: 0.01622572 - 0.019825899" optype="categorical" datatype="double"/>
    <DataField name="CODED1012" displayName="SPO2: 89.62968647 - 93.8" optype="categorical" datatype="double"/>
    <DataField name="CODED1013" displayName="DERIV1 SPO2: 0.002097087 - 0.015117552" optype="categorical" datatype="double"/>
    <DataField name="CODED1014" displayName="DERIV1 RESP RATE: -0.08696604 - -0.02060965" optype="categorical" datatype="double"/>
    <DataField name="CODED1015" displayName="DERIV1 RESP RATE CMOVAV: -0.00340129 - 0.001206744" optype="categorical" datatype="double"/>
    <DataField name="CODED1016" displayName="NN INTERVAL: 0.687847487 - 1.386134038" optype="categorical" datatype="double"/>
    <DataField name="CODED1017" displayName="DERIV2 NN INTERVAL CMOVAV: 0.014283031 - 0.046900061" optype="categorical" datatype="double"/>
    <DataField name="CODED1018" displayName="SPO2: 98.91534705 - 99.84375" optype="categorical" datatype="double"/>
    <DataField name="CODED1019" displayName="DERIV2 SPO2: 0.00061204 - 0.001579882" optype="categorical" datatype="double"/>
    <DataField name="CODED102" displayName="DERIV1 MULTI BICARB CMOVAV: 0.00008175 - 0.000082146" optype="categorical" datatype="double"/>
    <DataField name="CODED1020" displayName="DERIVT MULTI ION: -0.00109469 - -0.00002736" optype="categorical" datatype="double"/>
    <DataField name="CODED1021" displayName="DERIV1 RATIO CMOVAV: -0.00001006 - -3.4143E-51" optype="categorical" datatype="double"/>
    <DataField name="CODED1022" displayName="NN600: 1 - 3" optype="categorical" datatype="double"/>
    <DataField name="CODED1023" displayName="DERIV2 HR ECG 2: -1.00886608 - -0.00879864" optype="categorical" datatype="double"/>
    <DataField name="CODED1024" displayName="NN INTERVAL: 0.441642062 - 0.471441857" optype="categorical" datatype="double"/>
    <DataField name="CODED1025" displayName="DERIV2 MULTI ION CMOVSTD: 1.364081E-7 - 2.547968E-7" optype="categorical" datatype="double"/>
    <DataField name="CODED1026" displayName="DERIV2 SPO2 CMOVAV: 0.00077328 - 0.000186636" optype="categorical" datatype="double"/>
    <DataField name="CODED1027" displayName="DERIV2 PLATELET COUNT: -2.94728E-9 - -1.5087E-9" optype="categorical" datatype="double"/>
    <DataField name="CODED1028" displayName="PNN400: 0.02293578 - 0.142857143" optype="categorical" datatype="double"/>
    <DataField name="CODED1029" displayName="DERIV2 HEART RATE CMOVAV: -0.00016189 - -0.00004036" optype="categorical" datatype="double"/>
    <DataField name="CODED103" displayName="MULTI BASE CMOVSTD: 6.397148302 - 6.444158442" optype="categorical" datatype="double"/>
    <DataField name="CODED1030" displayName="DERIV1 SPO2: 0.015118904 - 0.037379534" optype="categorical" datatype="double"/>
    <DataField name="CODED1031" displayName="NUM SECS: 118 - 129" optype="categorical" datatype="double"/>
    <DataField name="CODED1032" displayName="PNN500: 0.3233709163 - 0.610859729" optype="categorical" datatype="double"/>
    <DataField name="CODED104" displayName="DERIV2 UREA NITROGEN CMOVAV: -2.6195E-10 - -1.2848E-10" optype="categorical" datatype="double"/>
    <DataField name="CODED105" displayName="MULTI LACTATE: 1.5 - 1.8" optype="categorical" datatype="double"/>
    <DataField name="CODED106" displayName="DERIV2 MULTI BASE CMOVAV: 5.727657E-08 - 6.723668E-8" optype="categorical" datatype="double"/>
    <DataField name="CODED107" displayName="DERIV2 MULTI-PH VEN: -1.1534E-10 - -2.980E-11" optype="categorical" datatype="double"/>
    <DataField name="CODED108" displayName="MULTI TON CMOVAV: 4.093846154 - 4.45624" optype="categorical" datatype="double"/>
    <DataField name="CODED109" displayName="MULTI TON: 4.69 - 4.77" optype="categorical" datatype="double"/>
    <DataField name="CODED11" displayName="DERIV1 MULTI PH CMOVAV: -9.68462E-7" optype="categorical" datatype="double"/>
    <DataField name="CODED110" displayName="DERIVT MULTI PH VEN CMOVSTD: 2.009864E-6 - 2.622579E-6" optype="categorical" datatype="double"/>
    <DataField name="CODED111" displayName="DERIV1 MAGNESIUM CMOVAV: 6.650743E-7 - 1.76845E-6" optype="categorical" datatype="double"/>
    <DataField name="CODED112" displayName="DERIV2 MULTI LACTATE CMOVAV: 1.33866E-10 - 6.68725E-10" optype="categorical" datatype="double"/>
    <DataField name="CODED113" displayName="DERIV1 MULTI LACTATE: -9.73923E-6 - -2.46708E-6" optype="categorical" datatype="double"/>
    <DataField name="CODED114" displayName="MULTI PO2 VEN CMOVSTD: 4.855041562 - 8.401710813" optype="categorical" datatype="double"/>
    <DataField name="CODED115" displayName="DERIV2 MAGNESIUM: 7.53241E-10 - 1.910706E-9" optype="categorical" datatype="double"/>
  </DataDictionary>
</PMML>

```

FIG. 18B

28/32

**FIG. 19**






FIG. 20



30/32

**FIG. 21**

31/32

**FIG. 22**

32/32

**FIG. 23****FIG. 24**