## (12) STANDARD PATENT

(11) Application No. AU 2012265587 B2

# (19) AUSTRALIAN PATENT OFFICE

(54) Title

#### PROSTHETIC FEMORAL STEM FOR USE IN HIGH IMPACT HIP REPLACEMENT

(51) International Patent Classification(s)

**A61F 2/32** (2006.01)

**A61B 17/74** (2006.01)

(21) Application No: **2012265587** (22) Date of Filing: **2012.12.17** 

(30) Priority Data

(31) Number (32) Date (33) Country 13/533,740 2012.06.26 US 13/366,156 2012.02.03 US 61/576,792 2011.12.16 US

(43) Publication Date: 2013.07.04
 (43) Publication Journal Date: 2013.07.04
 (44) Accepted Journal Date: 2015.02.12

(71) Applicant(s) Chow IP, LLC

(72) Inventor(s)

Chow, James C.; Rindal, Bjorn N.; Slater, Nicholas; Butters, Joshua

(74) Agent / Attorney

Fisher Adams Kelly Pty Ltd, L 29 12 Creek Street, Brisbane, QLD, 4000

(56) Related Art

US 2005/0251145 A1 (Desarzens et al.) 10 November 2005

## ABSTRACT OF THE DISCLOSURE

A total hip femoral prosthesis provides a shortened neck and more proximal coating level. Additionally, a means for strengthening the neck may be incorporated to increase stress transfer to the native bone. A method for reducing potential impingement that results from a high neck resection level includes instrumentation for cutting chamfered surfaces into the bone surrounding the proximal end of the implant. Instrumentation and methods of creating chamfered cuts are disclosed.




Fig. 7

AUSTRALIA	
P	Patents Act 1990

# ORIGINAL COMPLETE SPECIFICATION STANDARD PATENT

Invention Title:

"PROSTHETIC FEMORAL STEM FOR USE IN HIGH IMPACT HIP REPLACEMENT"

The following statement is a full description of this invention, including the best method of performing it known to me/us:

# PROSTHETIC FEMORAL STEM FOR USE IN HIGH IMPACT HIP REPLACEMENT

#### **BACKGROUND OF THE INVENTION**

This disclosure relates to prostheses and methods for total hip joint replacement. Specifically, this disclosure relates to prostheses that allow for high impact hip joint replacement and to instrumentation and methods for preparing the calcar region of the proximal femur that result in greater bone preservation. A high impact stem design is described that may result in greater bone preservation of the calcar region of the proximal femur, which results in better load sharing and greater overall capability. To ensure that the increase in bone preservation does not result in the limitation of the resulting range of motion, additional instrumentation and techniques are presented that allow the surgeon to create chamfer cuts on the calcar to minimize the potential for impingement.

Total hip replacement procedures seek to replace a hip joint that has deteriorated in its functionality, range of motion, weight bearing and most, if not all, other performance and lifestyle attributes. This degeneration is commonly found in patients that suffer from osteoarthritis, rheumatoid arthritis or osteonecrosis (avascular necrosis). Total hip replacement typically involves amputation of the femoral head, neck and a portion of the top of the femur in order to replace these structures with prosthetic components.

Individual skeletal development and postures vary from person to person. This is in part due to the three dimensional orientation of the hip socket relative to the proximal femur. The distance between the center of rotation of the femoral head and a reliable anatomical landmark, such as the lesser trochanter, may be described as the vertical offset of the head center from the lesser trochanter. This distance may be measured parallel to the femoral shaft axis, and is relevant to postoperative leg length. In the anterior-posterior view (AP view), the distance between the head center and the shaft axis may be described as the lateral offset of the shaft axis from the head center, or as the medial offset of the head center from the shaft axis. It is often referred to simply as "offset." Lateral offset may be relevant to postoperative hip abductor

function. Lateral offset is independent of the neck-shaft angle. However, lateral offset may be expressed in terms of the neck-shaft angle and the neck length, which is the distance along the neck axis between the head center and the shaft axis or some other reliable landmark. The neck-shaft angle varies through a range of angles, approximately 127-140 degrees for most people. The neck length varies as well. In the lateral view, the distance between the head center and the shaft axis may be described as the anteversion offset of the head center from the shaft axis, or simply anteversion, if the head center is anteriorly displaced from the shaft axis, or as retroversion if the head center is posterior to the shaft axis. Anteversion or retroversion may be relevant to postoperative range of motion.

Developed by Sir John Charnley in the 1960s, the original Charnley total hip arthoplasty was a stainless steel femoral prosthesis with a small collar, a rectangular cross section and a 22-mm femoral head. Subsequent designs evolved to include different head sizes (22, 25, 25.4 28, 32 and 35 mm), different femoral component lengths (ranging from 110 mm to 160 mm for standard prostheses), different cross sections (square, round, oval, I-beam) and a porous coating for bone ingrowth attachment and metal backing for the acetabulum (cemented or porous coated).

The neck-shaft angle and/or neck length of a prosthesis can also be highly varied in order to replicate natural anatomy or correct deformity. If the neck-shaft angle and/or neck length of a prosthesis are set so that the lateral offset is comparable to an average lateral offset value for intact normal femora, then this prosthesis may be said to have a conventional or standard offset. However, if the neck-shaft angle and/or neck length are set so that the lateral offset is relatively large, this is referred to as a high offset prosthesis. For example, the neck-shaft angle and/or neck length may be set so that the lateral offset is at least 10mm greater than the average lateral offset value.

Any of the relevant dimensions of a prosthesis may be set so that they are comparable to an average value for intact normal femora. These dimensions may be said to be conventional dimensions. Likewise, any of the relevant dimensions may be set incrementally greater than the corresponding average value, in which case the

prosthetic dimensions may be said to be augmented. It follows that smaller than average dimensions may also be selected.

A surgeon will typically measure both hip joints, including the neck-shaft angles, vertical offset, lateral offset, and leg length, prior to performing a total hip replacement procedure. These measurements allow the surgeon to match the replacement joint as closely as possible to the angles and dimensions of the original hip joint in order to achieve satisfactory range of motion, leg length, soft tissue tension, and stability. These measurements may also allow the surgeon to correct deformity or other conditions in and around the operated joint by matching the replacement joint to the contralateral hip joint.

#### STATEMENT OF INVENTION

In one aspect, the invention resides in a system for implanting a femoral hip prosthesis, comprising a broach, wherein the broach is shaped to be inserted into a proximal femur to form a hole, wherein the broach comprises a cutting guide contacting surface and a cutting guide connected to the broach, wherein the cutting guide comprises an engagement portion and a guide portion, wherein the engagement portion comprises a broach-contacting surface and a broach connection feature, wherein the broach-contacting surface contacts the cutting guide contacting surface, wherein the broach connection feature is complementary to the cutting guide connection feature, wherein the broach connection feature connects to the cutting guide connection feature; wherein the broach connection feature is on the broach-contacting surface, wherein the guide portion extends from the engagement portion, wherein the guide portion comprises a first guide surface, wherein the first guide surface defines a first resection plane.

In a further aspect, the invention resides in a system for implanting a femoral hip prosthesis, comprising a broach, wherein the broach is shaped to be inserted into a proximal femur to form a hole, wherein the broach comprises a cutting instrument contacting surface and a cutting instrument connection feature and a cutting instrument connected to the broach, wherein the cutting instrument comprises an engagement portion and a bone cutting surface, wherein the bone cutting surface

extends from the engagement portion, wherein the engagement portion comprises a broach-contacting surface and a broach connection feature, wherein the broach connection feature is complementary to the cutting instrument connection feature, wherein the broach connection feature connects to the cutting instrument connection feature, wherein the broach-contacting surface contacts the cutting instrument contacting surface, wherein the cutting instrument is rotatable, wherein when the instrument is rotated, the bone cutting surface engages and removes a section of bone.

In a yet further aspect, the invention resides in a system for total hip arthroplasty, the system comprising a hip prosthesis, the prosthesis comprising a stem, a neck, a means for strengthening the neck, and an exterior coating, wherein the stem comprises a longitudinal axis and a proximal portion, wherein the neck comprises a base and a head connecting feature, wherein the head connecting feature comprises a center point, wherein the means for strengthening the neck comprises additional material deposited on a superior surface of the stem, wherein the exterior coating comprises a proximal boundary, wherein a portion of the proximal boundary extends across a portion of the neck and an instrument for creating a chamfer cut in bone, wherein the instrument comprises a body, wherein the body comprises an engagement portion, wherein the engagement portion comprises a connecting feature, wherein the instrument comprises a cutting guide, wherein the cutting guide is coupled to the body.

In a still further aspect, the invention resides in a method for preparing a natural femur to receive a femoral hip prosthesis, comprising cutting across a neck of a natural femur to form a neck resection surface on a subcaptial portion of the neck, wherein the neck resection surface is bounded by an outer cortical edge, inserting a broach into a femoral metaphyseal cavity, wherein the outer cortical edge of the neck resection surface extends at least partially around the broach and chamfering the outer cortical edge of the neck resection surface.

Further aspects will be evident from the following description.

#### **BRIEF DESCRIPTION OF THE DRAWINGS**

Various examples of the present technology will now be discussed with reference to the appended drawings. It is appreciated that these drawings depict only typical examples of the invention and are therefore not to be considered limiting of its scope.

- FIG. 1 is a front view of a natural proximal femur, with dashed lines indicating anatomical variants with progressively larger lateral offsets;
- FIG. 2 is a front view of a prior art total hip femoral prosthesis, with dashed lines indicating a variant having a much larger lateral offset and a slightly larger vertical offset;
- FIG. 3 is a front view of another prior art total hip femoral prosthesis, with a cross-hatched area indicating a variant having a ninety degree neck angle, both prostheses having identical lateral and vertical offsets;
- FIG. 4 is a front view of yet another prior art total hip femoral prosthesis, with a cross-hatched area indicating a variant having the same lateral offset and a larger vertical offset;
- FIG. 5A is a front view of yet another prior art total hip femoral prosthesis, with a stem component, a conventional length neck, and a head component; and FIG. 5B is a front view of the stem and head of FIG. 5A with a longer neck component;
- FIG. 6A is a front view of a another prior art total hip femoral prosthesis with a short stem; and FIG. 6B is a front view of the femoral prosthesis of 6A inserted into a proximal femur with a traditional resection level;
- FIG. 7 is a cross section of a prosthetic hip stem with a shortened neck, high resection level inserted and high porous coating level into a femur;
- FIG. 8 is a cross section of a proximal femur with inserted hip stem with a traditional porous coating level;
- FIG. 9 is a side view of a proximal hip stem with a means for shortening the neck;
- FIG. 10A is a side view of a proximal hip stem with a means for shortening the neck; and FIG. 10B cross section of the superior surface of the proximal hip stem with a means for shortening the neck;
- FIG. 11A is side view of a hip stem with two different possible offsets; FIG. 11B is a table describing the geometric relationship between the diameter of the stem portion and the lateral offset;

FIG. 12A is a front view of a hip stem with a shortened neck, high neck resection level and a shortened stem; FIG 12B is a front view of the hip stem of 12A inserted into a proximal femur;

FIG. 13A is a cross section of a proximal femur that has been resected at the head neck junction; FIG. 13B is a top view of a proximal femur that has had the IM canal opened using a box chisel; FIG. 13C is a cross section of a femur that has been broached; FIG. 13D is a cross section of a femur with an inserted broach, before final neck resection; and FIG. 13E is a cross section of the proximal femur and broach of FIG. 13D after final neck resection;

FIG. 14A is a cross section of a proximal femur with a traditional resection level and a stem with a traditional porous coating level; FIG. 14B is a cross section of a proximal femur with a high resection level and a stem with a high porous coating level;

FIG. 15A is a superior view of a chamfered bone section with a high resection level; FIG. 15B is a front view of a chamfered bone section with a high resection level; FIG. 15C is a superior front view of a chamfered bone section with a lower resection level; FIG. 15D is a front view of a chamfered bone section with a lower resection level;

FIG. 16A is an isometric view of a proximal femur with a high resection level that has been broached; FIG. 16B is an isometric view of the proximal femur of FIG. 16A, with a shaded area indicating bone to be removed to form a chamfered anterior edge; FIG. 16C is a top view of the proximal femur with a high resection level and a chamfered anterior edge; FIG. 16D is a proximal femur with a high resection level and a chamfered anterior and posterior edge; FIG. 16E is a top view of a proximal femur with a high resection level and chamfered anterior edge, with a hip stem and a head; FIG. 16F is a top view of a rotated proximal femur with a high resection level and chamfered anterior edge, with a hip stem and acetabular cup inserted;

FIG. 17A is a side view of a broach that includes a mating feature for guiding a reamer along a trajectory to create a chamfer on the bone; FIG. 17B is a side

view of a planar reamer and a broach; FIG. 17C is a side view of a planar reamer engaged with a broach;

FIG. 18A is a bottom perspective view of a chamfered reamer with a cupshaped cutting feature; FIG. 18B is a cross section of the chamfered reamer of FIG. 18A;

FIG. 19A is a front view of a proximal femur with an inserted broach and protruding post; FIG. 19B is a front view of a chamfered reamer engaged with the broach of FIG. 19A;

FIG. 20A is a front view of a proximal femur with an inserted broach and protruding post after the calcar region has been reamed using a chamfered reamer; FIG. 20B is a front view of a proximal femur that has been reamed using a chamfered reamer and a hip stem has been inserted;

FIG. 21A is a superior view of a broach with a beveled edge; FIG. 21B is a top view of a broach with a beveled edge inserted into a proximal femur cutout with a beveled anterior edge;

FIG. 22A is a front view of a femur with a hip stem with a modular head mated with an acetabular cup and liner in a pelvis; FIG. 22B is a front view of a hip stem with a modular head inserted into an acetabular liner, where there is impingement between the bone and the acetabular implant;

FIG. 23A is a top perspective view of a broach inserted into a proximal femur; FIG. 23B is a top perspective view of a chamfer cutting guide attached to a broach inserted into a proximal femur;

FIG. 24 is a top perspective view of a chamfer cutting guide;

FIG. 25 is a bottom perspective view of a chamfer cutting guide;

FIG. 26 is a side view of a chamfer cutting guide;

FIG. 27A is a top view of a chamfer cutting guide attached to a broach inserted into a femur and a chamfered bone surface; and FIG. 27B is a top perspective view of a chamfer cutting guide attached to a broach inserted into a femur and a chamfered bone surface; and

FIG. 28 is a diagram illustrating methods of preparing anterior and/or posterior chamfer cuts.

#### **DETAILED DESCRIPTION**

The present disclosure sets forth a prosthetic high impact femoral hip stem including shortened neck and a means for strengthening the neck. Additionally, the present disclosure sets for a means for reducing potential impingement with the prosthetic hip stem, and maintaining the patient's natural range of motion.

In this specification, standard medical directional terms are employed with the ordinary and customary meanings. Superior means toward the head. Inferior means away from the head. Anterior means toward the front. Posterior means toward the back. Medial means toward the midline, of plane of bilateral symmetry, of the body. Lateral means away from the midline of the body. Proximal means toward the trunk of the body. Distal means away from the trunk. Additionally, standard hip anatomical terms are employed with their ordinary and customary meanings.

Referring to FIG. 1, an intact natural femur 110 is shown in solid lines in a front view, which may also be described as an anterior-posterior (AP) view. The femur 110 has a proximal end 100 which includes a shaft 102, a neck 104 and a head 106. The shaft 102 has a center longitudinal axis 112. The neck 104 has a center longitudinal axis 124 which is oblique to the shaft axis. The head 106 may be approximately spherical. The head has a center point 114, or center of rotation center 114 is spaced apart from the shaft axis 112 by a lateral offset distance 130, which is an average lateral offset value, otherwise known as a conventional offset or standard offset.

Lateral offset 130 may be approximately 39.8mm to 54.2mm, based on the work of Rubin, et al. as published in 1992 in the British Journal of Bone and Joint Surgery (JBJS BR), volume 74B, pages 28-32, which is incorporated by reference herein in its entirety. Alternately, lateral offset 130 may be 36.2mm to 49.8mm, based on the work of Noble, et al. as published in 1988 in Clinical Orthopedics and Related Research (CORR), Number 235, pages 148-165, which is incorporated by reference herein in its entirety.

FIG. 1 also illustrates an anatomical variant of femur 110, which has another approximately spherical femoral head 116 with a center point 118. The center 118, which may lie on the neck axis 124, is spaced apart from the shaft axis 112 by a

second lateral offset distance 132 which is greater than lateral offset 130. This distance has been exaggerated in FIG. 1 for clarity. This variant may be said to have a high lateral offset. The center 118 may also have a greater vertical offset than center 114 from any reliable anatomical landmark on femur 110, and may thus be offset from the neck axis 124. It can be appreciated that the neck 104 in this variant is naturally long enough to reach from the shaft 102 to the head 116.

FIG. 1 also illustrates another anatomical variant of femur 110, which has yet another approximately spherical femoral head 120 with a center 122. The center 122, which is shown lying on neck axis 124, is spaced apart from the shaft axis 112 by a third lateral offset distance 134 which is greater than the second lateral offset 132. This distance has also been exaggerated in FIG. 1 for clarity. This variant may be said to have a high offset. More specifically, the offset may be characterized as an extrahigh offset. The center 122 may also have a greater vertical offset than center 118 from any reliable anatomical landmark on femur 110. It can be appreciated that the neck 104 in this variant is naturally long enough to reach from the shaft 102 to the head 120.

Rubin (JBJS BR 1992) reports offsets up to 62.8mm. Offsets up to 68.6mm may be predicted from Rubin's mean of 47.0mm plus three standard deviations of 7.2mm. Noble (CORR 1988) reports offsets up to 61.0mm. Offsets up to 63.4mm may be predicted from Noble's mean of 43.0mm plus three standard deviations of 6.8mm.

Referring to FIG. 2, a prior art total hip femoral prosthesis 200 is shown implanted in the proximal end 100 of femur 110. Modular femoral prosthesis 200 includes a stem component 206, a neck component 208, and a head component 210. The neck 208 and head 210 and stem 206 are secured together. The stem 206 is shown within a prepared canal in the femoral shaft 102. A femoral neck osteotomy has been performed, removing the neck 104, head 106, and a portion of the proximal end 100, leaving an angled resection surface 108. The stem 206 has a center longitudinal axis 214 which may align with the shaft axis 112 when the stem is implanted in the shaft 102. The neck 208 has a center longitudinal axis 216 which may be oblique to the stem axis 214. The head 210 is spherical and has a center point 212 which lies on the

neck axis 216. The center 212 is spaced apart from the stem axis 214 by a lateral offset distance 240, which is a conventional lateral offset.

FIG. 2 also illustrates another prior art total hip femoral prosthesis 220, which includes the stem 206, another neck component 222, and another head component 224. The neck 222 has a center longitudinal axis 230 which may be oblique to the stem axis 214. The head 224 is spherical and has a center point 226 which lies on the neck axis 230. The center 226 is spaced apart from the stem axis 214 by a lateral offset distance 242 which is greater than lateral offset 240. Lateral offset 242 may be called a high offset compared to lateral offset 240. The center 226 may also have a greater vertical offset than center 212 from any reliable anatomical landmark on femur 110. Neck 222 is approximately 38mm long measured from base to end along the neck axis 230.

It can be appreciated that high offset neck 222 is longer than conventional neck 208. Therefore, neck 222 may experience higher stresses than neck 208 during use. Neck 222 may bend, crack, or break due to high service stresses and/or accidental overload from trips, slips, falls, or other trauma, particularly over years *in vivo*. Neck 222 may be expected to fail in or near the indicated region 228. Failure of the neck 222 will require corrective surgery to replace at least the broken neck, with attendant rehabilitation and risks of significant complications.

Referring to FIG. 3, a prior art total hip femoral prosthesis 302 is shown implanted in the proximal end 100 of femur 110. Modular femoral prosthesis 302 includes a stem component 304, a neck component 306, and a head component 308. The neck 306 and head 308 and stem 304 are secured together. The stem 304 is shown within a prepared canal in the femoral shaft 102. The femoral shaft 102 has been prepared as set forth above. The stem 304 has a center longitudinal axis 310 which may align with the shaft axis 112 when the stem is implanted in the shaft 102. The neck 306 has a center longitudinal axis 314 which may be oblique to the stem axis 310. The head 308 is spherical and has a center point 312. The center 312 is spaced apart from the stem axis 310 by a lateral offset distance 340.

FIG. 3 also illustrates another prior art total hip femoral prosthesis 320, which includes a longer stem component 322, a horizontal neck component 324 and

another head component 326. The stem 322 is aligned with stem axis 310. The neck 324 has a center longitudinal axis 328 which may be perpendicular to the stem axis 310, or nearly so. The head 326 is in the same position as head 308. Therefore, head 326 has the same lateral offset 340 and vertical offset as head 308.

It can be appreciated that neck 324 may be shorter than neck 306, however, stem 322 and neck 324 protrude proximally from femur 110 more than stem 304 and neck 306 do. Therefore, stem 322 and neck 324 may hinder abduction by impinging surrounding anatomical structures, such as the lateral pelvis.

Referring to FIG. 4, a prior art total hip femoral prosthesis 402 is shown implanted in the proximal end 100 of femur 110. Modular femoral prosthesis 402 includes a stem component 403, a neck component 404, and a head component 406. The neck 404 and head 406 and stem 403 are secured together. The stem 403 is shown within a prepared canal in the femoral shaft 102. The femoral shaft 102 has been prepared as set forth above. The stem 403 has a center longitudinal axis 410 which may align with the shaft axis 112 when the stem is implanted in the shaft 102. The stem 403 has a placement slot 408 for a stem inserter instrument (not shown), a tool used for driving the stem into the prepared canal within the femur. The neck 404 has a center longitudinal axis 414 which may be oblique to the stem axis 410. The head 406 is spherical and has a center point 412. The center 412 is spaced apart from the stem axis 410 by a lateral offset distance 440.

FIG. 4 also illustrates another prior art total hip femoral prosthesis 422, which includes another stem component 424, another neck component 426, and another head component 428. The stem 424 is aligned with stem axis 410. The neck 426 has a center longitudinal axis 432 which may be oblique to the stem axis 410. The head 428 is spherical and has a center point 430. The center 430 is spaced apart from the stem axis 410 by lateral offset distance 440. It can be appreciated that center 430 is farther than center 412 from any landmark on femur 110 in a direction parallel to axis 410. Prosthesis 422 may be said to have a high vertical offset compared to prosthesis 402.

Referring to FIG. 5A, a prior art modular total hip femoral prosthesis 500 includes a stem component 502, a neck component 506, and a head component 508.

The neck 506 and head 508 and stem 502 are secured together. The stem 502 has a center longitudinal axis 512 which may align with the shaft axis 112 when the stem is implanted in the shaft 102. The stem 502 also has a pocket 504 which receives a distal portion of the neck 506 when the prosthesis 500 is assembled. The neck 506 has a center longitudinal axis 516 which may be oblique to the stem axis 512. A distal center of rotation 510 may be defined on the neck axis 516 in the distal portion of the neck 506. The distal center of rotation 510 may be in the pocket 504 when the neck 506 is secured to the stem 502. The head 508 is spherical and has a center point 514. The center 514 may coincide with the neck axis 516 in a proximal portion of the neck 506 when the head 508 is secured to the neck 506, thus the center 514 may also be referred to as a proximal center of rotation of the neck. The center 514 is spaced apart from the distal center 510 by a medial-lateral distance 518 and a distance 520 along the neck axis 516. Medial-lateral distance 518 is approximately 20mm to 50mm. Distance 520 is approximately 20mm to 30mm.

Referring to FIG. 5B, a prior art high offset modular total hip femoral prosthesis includes stem 502 and head 508 with a longer modular neck component 550. The neck 550 has a center longitudinal axis 552 which may be oblique to the stem axis 512. The center of head 508 in this arrangement is at point 554. The center 554 is spaced apart from the distal center 510 by a medial-lateral distance 522 and a distance 524 along the neck axis 552. Distance 524 is greater than distance 520. Medial-lateral distance 522 is greater than medial-lateral distance 518 by the increment 526. Increment 526 is approximately 10mm to 20mm. Therefore, medial-lateral distance 522 is approximately 30mm to 70mm. Distance 524 is approximately 30mm to 40mm. Disadvantageously, the increase in lateral offset places significantly greater stress and strain on the longer modular neck 550, as described for neck 222.

Referring to FIGS. 6A and 6B, a standard, prior art short wedge-style prosthetic hip stem is illustrated with a traditional neck resection level. As illustrated in FIG. 6A, the short wedge-style hip stem 800 includes a traditional porous coating level 802 and a traditional effective neck length 804, as defined by the distance from a proximal boundary 806 of the porous coating level to a center point 807 of a head connection feature 808. The short wedge-style prosthetic hip stem 800 includes a

shortened stem portion 810. The shortened stem portion 810 includes an elbow feature 812, which may also be described as a trochanteric relief transition. A distal boundary 805 of the porous coating level is shown below the trochanteric relief transition 812. Referring to FIG. 6B, the standard short wedge-style hip stem 800 is shown inserted into a proximal femur with a traditional neck resection level 814. The neck resection is level corresponds with the proximal boundary 806 of the porous coating.

Having described various commercially available styles of total hip prosthesis, the following description relates to examples of a total hip prosthesis that includes a shortened neck and more proximal coating allowing for a higher neck resection. Another example includes a means for strengthening the neck. Additionally, the following description provides a method for preparing the calcar region such that potential impingement is reduced, allowing the patient to maintain a normal range of motion.

Referring to FIG 7, an example of a prosthetic femoral hip stem 600 with a shortened neck and more proximal coating is shown implanted into a proximal femur. Femoral hip stem 600 may include features similar to those described in FIGS. 1-6, including a conventional lateral offset 601, a standard stem length and a standard proximal body design. The stem component 602 is shown to extend into the intramedullary canal. Femoral hip stem includes a head connection feature 606 shaped to receive a modular femoral head. Head connection feature 606 may include a center point 605. Additionally, femoral hip stem 600 may be designed for a higher neck resection than is used with conventional stems, allowing for a shortened prosthetic neck 604. The neck 604 may be medialized in this example. The neck axis has been shifted away from the stem axis, compared to a conventional stem. Portions of the stem 602, neck 604 and head connection feature 606 may be fabricated from different materials according to the requirements of each portion.

Femoral hip stem 600 may also include a superior surface 620 that extends along a transverse plane that is proximal to the stem 602, and may intersect a lateral portion of stem 602 and neck 604. In other words, superior surface 620 may blend the superior portions of the stem 602 and a base portion of the neck 604 together. The

superior surface may by smooth and planar, or may include surface texture elements such as roughening, and may also be contoured.

The shortened neck 604 may be fabricated with a substrate of one material, and include an exterior coating 608 of another material. The coating 608 may be porous to enhance bone ingrowth around the neck 604. The shortened neck 604 may allow for a high level of porous coating 608, wherein at least a portion of a proximal boundary 610 of the more proximal coating 608 may extend from a proximal portion of the stem 602 along a plane parallel to the superior surface 620. Bone ingrowth onto the shortened neck may facilitate load sharing with the shortened neck, allowing for a high impact application. The proximal boundary 610 of the more proximal coating 608 may then extend distally at an angle across a base portion of the neck 604. The angled portion of the proximal boundary 610 may be located at a first distance 612 from the center point 605 of the head connecting feature 606.

The lateral extent of shortened neck 604 may be defined as the lateral distance between the resection level on the stem, or the porous coating level on the stem, and the center of the head connecting feature 606, or the center of an actual head, if present. More specifically, the lateral extent may be defined as the lateral distance between the intersection of the neck axis and the resection or porous coating level, and the center of the head connecting feature 606, or head. The lateral extent of the shortened neck 604 may be a percentage of the total lateral offset distance described previously, for example with reference to FIG. 1. The lateral extent may be about 50% of the lateral offset, or it may be smaller, such as 40%, 30% or 20% of the total medial offset. For example, if the total lateral offset is 45 mm, the lateral extent may be defined as 30% of the total lateral offset, or 15 mm. It can be appreciated that a smaller lateral extent correlates with a shorter neck, as measured along the neck axis.

Comparatively, a femoral hip stem that is similar to hip stem 600, but with a longer neck and a traditional porous coating level is shown in FIG.8. The traditional porous coating has a proximal boundary is located a second distance 614 from the center point 605 of head connecting feature 606, wherein the second distance 614 is larger than the first distance 612 shown in FIG. 6.

By shortening the neck 604 and including a more proximal coating 608, as shown in FIG. 7, the amount of unsupported neck length may be reduced and the moment arm length may be decreased, as a more proximal coating may enhance bony ingrowth closer to the head portion 606, thus reducing the stress put on the neck 606 by increasing the amount of load that can be shared with native bone. The high impact hip stem described in FIGS. 7 and 9-11 is optimized for larger diameter heads (such as ≥36 mm), providing a range of motion comparable to the natural hip joint. By designing the prosthetic neck geometry around a larger diameter prosthetic head, the prosthetic neck length can be shortened while maintaining a desirable ratio between the diameter of the prosthetic head and the diameter of the neck at the level where acetabular impingement would occur.

Another example of a high impact femoral hip stem 600 may include a means for strengthening the neck. Referring to FIGS. 9 and 10A-10B, a proximal portion of femoral hip stem 600 is shown. It can be appreciated that stem 600 may include additional material deposited on the superior surface 620, relative to earlier designs, to create a rounded superior surface 622 that serves as a means for strengthening the neck. The addition of material to the superior surface to create a means for strengthening the neck may include material that is the same as the material used to fabricate the neck and stem portions, or may be a different material. FIGS. 10A and 10B provides a cross section A-A of the rounded superior surface 622, where additional material is present on the superior surface 620. Alternatively, the additional material may create a proximally elevated cross-section that is polygonal, rectangular or otherwise irregularly shaped.

The rounded superior surface 622 may serve to maintain a constant distance between the intended femoral head center location and the intersection of the lateral neck and the stem superior surface. This "effective neck length" 624 may remain approximately constant as the size of the overall implant is varied (stem diameter and offset). The additional material 620 may thus provide a constant, effective neck length 624, which is defined as the length between the intersection of the neck 604 and the additional material 620 and a center point of the head component 606. The porous coating 608 may also extend over the rounded superior surface 622.

The designs herein, such as that illustrated in FIGS. 9-10B, may be well adapted for prosthetic femoral heads larger than 22 mm, preferably larger than 32 mm, and preferably larger than 36 mm. Larger prosthetic heads may provide relatively greater range of motion than do smaller heads. The increased range of motion may be enjoyed on its own merits, or it may influence other choices in the overall design of a large-head total hip arthroplasty system. For example, instead of providing an actual increase in range of motion in the large-head system, the neck diameter or width may be augmented to improve its strength, a higher neck resection may be used, and/or a thinner acetabular component may be provided. A higher neck resection may permit the design to incorporate a higher level of porous coating, a shorter neck, and/or a medialized neck. While these choices may result in a large-head system whose range of motion is comparable to that of a small-head system, they may also result in a system whose strength and load transfer characteristics are better adapted for high demand patients.

Referring to FIGS. 11A and 11B, an additional example of the high impact femoral hip stem may combine the shortened neck length 604 described in FIG. 7 with the rounded superior surface 622 described in FIGS. 8 and 9, as a means for strengthening the neck. By combining these features, the neck may be strengthened enough such that the maximum femoral offset available for a given stem diameter or size may be increased. Specifically, this may be advantageous for patients with small intramedullary canal anatomy and large lateral offsets. This femoral morphology may be referred to as a "champagne flute" morphology. Referring to FIG. 11A, two different possible offsets are illustrated- a high offset 642, which may have lateral offsets comparable to the high offset offerings of standard stem designs, and a high offset plus 640 design, which may have increased lateral offsets.

The table presented in FIG. 11B presents an example of a geometric relationship between the diameter of the stem portion 602 and the lateral offset. The leftmost column of the table lists various distal stem diameters. The rightmost two columns present lateral offset values for the two designs. The rightmost column illustrates a 5 mm increase in lateral offset, although other increments are contemplated. For example, the high offset plus design may incorporate lateral offsets

which are greater than one or two standard deviations above an average lateral offset for a patient population. Lateral offsets which are greater than three or more standard deviations above the average are also contemplated.

Referring to FIGS. 12A-12B, an alternative embodiment of a high impact femoral hip stem with a shortened stem is illustrated. Femoral hip stem 850 may be similar to the hip stem described in FIGS. 6A and 6B, but may include a more proximal porous coating 852, a higher neck resection level 854, and thus, a shorter "effective neck length". Femoral hip stem 850 may include a shortened stem 856 with a trochanteric relief transition 858, and a medial curvature 860 with a continuous flare and gradual curvature. Femoral hip stem 850 may also include a head connecting feature 862, similar to previous embodiments described herein.

The shortened stem 856 may be less than 6 cm in length, or may be longer. The more proximal porous coating 852 may include a proximal boundary 851 and a distal boundary 853. The proximal boundary 851 of the more proximal porous coating 852 is illustrated to lie closer to the head connecting feature 862 than the proximal boundary 806 of the porous coating illustrated in FIGS. 6A and 6B. Also, the distal boundary 853 of the more proximal porous 852 coating may be located at or proximal to the trochanteric relief transition 858. By placing the distal boundary 852 of the more proximal porous coating 852 at this location, the implant may be more easily removed during a revision surgery, as the coating does not extend around the bend 858.

Further, when femoral hip stem 850 is inserted into a proximal femur, as illustrated in FIG. 12B, the neck resection level 854 is higher than the neck resection 814 shown in FIG. 6B. By shifting the neck resection level proximally, and thus shifting the porous coating proximally, subsequent revisions to the hip stem may become easier when compared to traditional hip stems with traditional porous coating levels. Also, similar to previous embodiments, the high neck resection level may allow for greater preservation of the natural femur, specifically in the calcar region, allowing the stem to share more of the load with the native femur and increasing the overall strength of the construct. Bony integration may be achieved closer to the prosthetic head. Because the porous coating is more proximal, and higher up on the medial flare, a greater compressive stress may occur in the bone at the interface, which may help to

facilitate greater bone ingrowth. This may further reduce the unsupported length of neck. Overall, this may substantially prevent loosening of the stem in the femur, reduce migration of the stem within the femur, and provide greater implant stability.

Femoral hip stem 850 may also include a means for strengthening the neck, as described in previous embodiments.

Any of the stem components set forth herein may be fabricated from biocompatible materials. Examples of biocompatible materials include, but are not limited to, metals such as stainless steels, titanium and its alloys, cobalt-chromemolybdenum alloys, and other chromium alloys; polymers such as polyetheretherketone (PEEK) and polyaryletherketone (PAEK); ceramics such as alumina and zirconia; and composites such as carbon fiber reinforced epoxy resin. Portions of a stem component may be fabricated from different materials according to the requirements of each portion. A stem component may be fabricated with a distal portion of one material and a proximal portion of another material. A stem component may be fabricated with a substrate of one material and a coating of another material. The material or materials may be solid (i.e., non-porous) or porous. The surface of a stem component may be smooth or rough on a micro- or macroscopic level. A stem component may be fabricated by casting, forging, machining, or a combination of methods. Castings may be applied by thermal, chemical, electrical, or comparable means. For example, coatings may be sintered, sputtered, vapor deposited, ion implanted, electroplated, and the like.

Any of the neck components set forth herein may be fabricated as described above for the stem components. Neck components may be preferably fabricated with materials, surface treatments, and methods of fabrication which enhance strength, or minimize unavoidable reductions in strength.

FIGS. 13A-13E illustrate steps in a method for preparing and inserting the high impact hip stem with a high final resection level. An osteotomy may be performed by removing the head to create a prepared resection surface at or near the head-neck junction as illustrated in FIG 13A. The intramedullary (IM) canal may then be opened, using a box chisel or other comparable surgical instrument as shown in FIG. 13B.

Once the IM canal has been opened, as shown in FIG. 13C the surgeon may then ream the distal IM canal and broach the proximal IM canal and proximal metaphysis, as illustrated in FIGS. 13D and 13E. Next, the neck resection may be completed using the broach as the resection guide for a planar calcar reamer or saw, FIG. 13E.

FIGS. 14A and 14B provide a comparison between a traditional resection level (FIG. 124) and the final resection level according to the technique described in FIGS. 11A-11E. The final resection level (FIG. 14B) is higher on the femoral neck than traditional resectioning techniques, resulting in greater load sharing with native bone.

The high impact stem design described in FIGS. 7 and 9-11 may use the surgical technique described in FIGS. 13A-13E for greater bone preservation of the calcar region of the proximal femur during total hip arthroplasty, and thus enhanced load sharing. However, the increase in bone preservation potentially limits the patient's range of motion due to impingement.

Previous designs utilize a lower femoral neck cut and alter the implant neck geometry, including the incorporation of neck flats, to provide relief against acetabular component impingement when using smaller head sizes. When using larger head sizes, this same geometry alteration can be made in the surrounding bone, for example, by using chamfer cuts, relief, rounding, edge break, or other modifications in the zone of potential impingement. Chamfer cuts or other cuts made in the surrounding bone may provide relief where impingement will occur against the acetabular component. This may allow for a thicker implant neck geometry to be used, since impingement relief is addressed in the bone and not the implant. Any method can be used to alter the bony geometry to this effect, including removing portions of bone to reveal chamfered, planar surfaces. Some proximal bone may be retained for load sharing in a high impact application, such as a bony "shelf" to help transfer load from the neck portion of an inserted hip stem. Alternative bone surface geometries may include a rounded surface, an L-shaped cut-out or otherwise contoured surface resection geometry.

Figures 15A and 15B illustrate an example of a chamfered bone section that may be used in conjunction with the high impact hip stem that has previously

been described. The chamfered bone surfaces shown in FIGS. 15A and 15B may maintain a high level of calcar bone, while allowing for a typical range of motion. Figures 15C and 15D illustrate another example of a chamfered bone, in which the resection surface is more angulated, while still retaining a high level of calcar bone, and also allow for a typical range of motion.

FIGS. 16A-25B illustrate instruments and methods that may be used for creating the chamfered surfaces to reduce impingement and allowing for a hip stem with a strengthened neck. FIGS. 16A-16F illustrate additional steps to further prepare the anterior calcar chamber to increase range of motion and to prevent possible impingement of the prosthetic components and/or host bone. FIG. 16A illustrates a proximal femur that may have been prepared according to the technique of FIGS. 13A-13E. A portion of the femoral head has been removed to reveal a horse-shoe shaped cavity 650 that may allow for access to the IM canal. The depth and width of the cavity may be variable, according to the patient's anatomy. Alternatively, the cavity may be rectangular or irregularly shaped. In FIG. 16B, the proximal femur has been further prepared by introducing a beveled or chamfered anterior edge 652 to the cavity 650. FIG. 16C provides a top view of the cavity 650, which may include a first, non-beveled posterior edge 654 and a beveled anterior edge 652. Alternatively, both the anterior and posterior edge may be beveled. In this case, the posterior and inferior bone stock may be preserved for maximal stress transfer, as illustrated in FIG. 16D. When the anterior edge is beveled, the range of motion may be enhanced, as shown in FIGS. 16E and 16F. Referring to FIGS. 16E-16F, a top view of a proximal hip stem inserted into a femur through the cavity 650 is shown, with the head connecting component 606 connected to a modular head 665, which is received by an acetabular cup 666. The implanted hip stem may be the hip stem 600 described in previous examples. FIG. 16E illustrates the hip stem 600 and connected modular head 665 in a neutral, non-rotated position. FIG. 16F illustrates the hip stem 600 and connected modular head 665 rotated such that the anterior, beveled edge 452 of the cavity 450 comes into contact with the acetabular cup 666. It can be seen that the beveled edge 652 reduces the potential for impingement by allowing the proximal femur and to rotate further before contacting the anterior edge of the acetabular cup 666.

Referring to FIGS. 17A-17C, a beveled anterior edge 652 may be prepared by using a planar reamer 670, sometimes referred to as a calcar planer or calcar reamer, on a post 672 engaged with a broach 674. FIG. 17A shows a broach 674. The broach includes an oblique proximal hole 676 which is perpendicular to the desired angle of the beveled anterior edge 652. With the broach 674 fitted into the proximal femur, the post 672 is inserted into the hole 676 to align the reamer 670 with the hole. As the reamer 670 rotates, as indicated by motion arrows 677 and advances, it cuts away the proximal anterior edge to form the bevel 652. The reamer 670 may include a central cavity sized and shaped to receive a proximal anterior corner or shoulder of the broach without interference while the reamer is spinning.

The broach may include only one hole, in which case separate right and left broaches may be provided. The broach may also include two holes so that a single broach is suitable for right and left femurs.

Referring to FIGS 18A-20B illustrate another example of a method and instrument that may be used to prepare a femoral hip stem for a hip prosthesis with a high resection level, while retaining a normal range of motion in a patient. Referring to FIGS. 18A and 18B, a chamfered reamer is shown. Chamfered reamer 900 includes an elongated handle portion 902, connected to a cup-shaped cutting surface 904. The cupshaped cutting surface may also be referred to as a bone cutting surface. The cupshaped cutting surface 904 may be conical, or include varying degrees of curvature. The elongated handle portion 902 is connected to the cup-shaped cutting portion 904 at a shoulder 906. The shoulder 906 may have a greater diameter than the elongated handle 902. The cup-shaped cutting portion 904 may be located distal to the shoulder 906, and may flare outwards. The cup-shaped cutting portion 904 may include an inner surface 903, and an outer surface opposite 905 of the inner surface 903. The cupshaped cutting portion 904 may include a sharpened edge portion 908 at the intersection of the inner 903 and outer surfaces 905. The cup-shaped cutting portion 904 may include a plurality of apertures 910, which may also be referred to as windows, which may include sharpened edges 911 that facilitate removal of bone. The window edges 911 may otherwise be flattened or rounded. At the apex of the cupshaped cutting portion 904, the chamfered reamer 900 may include a female feature 912 shaped to receive a portion of a post that is connected to a broach. The female feature 912 may include a socket, or hole, that extends from the inner surface 903 of the cup-shaped cutting portion 904 into the shoulder 906. The mating feature 912 may otherwise be a complementary male feature shaped to be received by a groove or socket in a broach.

Referring to FIG. 19A, a broach with a post 920 is shown inserted into a proximal femur, after the head of the proximal femur has been resected at the headneck junction. The broach may have been inserted into the femoral cavity until the broach is flush with the resection level, and the post is oriented such that it extends substantially normal to the resection plane. The chamfered reamer 900 described in FIGS. 18A and 18B is shown mating with the inserted broach in FIG. 19B. When the chamfered reamer 900 is engaged with the broach post 920, the sharpened edge portion 908 may lie flush with the resection plane. The instrument 900 may then be actuated to ream the calcar section by rotating and applying force to the chamfered reamer to remove bone tissue.

Referring to FIG. 20A, a proximal femur is shown after it has been reamed using the chamfered reamer 900, with the broach 920 still implanted in the proximal femur. After reaming, the broach may be removed and a femoral hip prosthesis, such as those previously described, may be inserted, as illustrated in FIG. 20B.

Another example of a calcar reamer may include a central aperture shaped to receive a proximal portion of a broach that has been inserted into a proximal femur. The structure and method of use of the calcar reamer may be similar to the reamers shown in FIGS. 17A-20C and, although the calcar reamer may include at least one chamfered cutting surface, such that as the reamer is rotated, the chamfered cutting edge acts in a saw-like capacity to create a complimentary chamfered bone surface on one or both sides of the broach. These instruments may include a variety of cutting edge geometries to create differently shaped bone resections. Alternatively, a resection may be done free-handed to remove bone portions surrounding the hip stem and reduce potential infringement.

FIGS. 21A-21B illustrate another arrangement for preparing the beveled anterior edge 652 by using a saw guided by a beveled edge 678 of a broach 674. In

this example, the proximal anterior corner of the broach carries a beveled portion 678 which establishes the desired resection plane for the beveled anterior edge 652. With the broach 676 fitted into the proximal femur, the saw is aligned with and moved along the beveled portion 678. As the actuated saw advances, it cuts away the proximal anterior edge to form the bevel 652.

The broach may include one or two beveled portions 678.

FIGS. 22-27 illustrate another example of an instrument and method that may be used for preparing the proximal femur during total hip arthroplasty to minimize the potential for impingement and maintain the patient's natural range of motion. FIGS. 22A and 22B illustrate a prosthetic femoral hip stem connected to a modular head inserted into a femur. The hip stem may be the high impact hip stem 600 described previously in this application. Referring to FIGS. 22A and 22B, no chamfer resection cuts have been performed after the hip stem 600 was inserted into the femur. The connected modular head portion can be seen to interface with an acetabular component. It can be seen in FIG 22B that without chamfer resection cuts, there may be posterior impingement 681 of the joint at the bone-acetabular cup interface, thus limiting the patient's range of motion.

FIGS. 23-27 illustrate a chamfer cutting guide 700 that may be fitted or otherwise attached to the proximal portion of a broach 674 that has been inserted into the proximal femur, and may be used to create chamfered edges in the native femur that reduce impingement of the prosthetic head against the anterior and posterior portions of the natural femur.

Referring to FIG. 23A, a broach 674 is shown inserted into the proximal femur. The broach 674 may be inserted into the proximal femur using a method similar or identical to the process described previously in this application. FIG. 23B illustrates a chamfer-cutting guide 700 that may be attached to a proximal portion of the broach 674.

FIG. 24 provides a perspective view of chamfer cutting guide 700, which may include a first portion 710, which may be referred to as an engagement portion and a second portion 750, which may be referred to as a guide portion. The engagement portion 710 may act to secure chamfer cutting guide 700 to a broach 674

or other anchor feature in bone, and may include a proximal surface 712, and a distal broach contacting portion 714 opposite the proximal surface. The proximal surface 712 may extend between a rounded edge 716 and a guide portion interface 718.

The distal broach contacting 714 portion, as illustrated in FIG. 25, may have a step-like configuration, and include a first step 720, a second step 722 that is offset distally from the first step 720 and a third step 724 that is offset distally from the second step 722. The number and orientation of steps may be variable. The distal broach contacting surface 714 may also include various alternative configurations shaped to interface with the bone portion surrounding the inserted broach, specifically, the bone-contacting portion may interface with the bone on the interior portion of the cavity 650 described in FIGS. 16A and 16B. Alternative configurations may include a flat, sloped, or otherwise contoured surface to match the surface geometry of the patient's bone.

The surface of the first step 720 on the distal broach contacting portion 714 may include at least one aperture 726, and in this case, includes two relatively small and substantially circular apertures.

The first step 472 and second step 724 may be connected by a first vertical surface 728, which may intersect the first step 720 at a rounded concave corner 730. The corner may otherwise by sharp, beveled, coved or eased. The second step 724 may include at least one larger aperture 732, as shown best in FIG. 25. The larger aperture 732 shaped to receive a proximal portion of a broach. The aperture 732 may extend entirely through the engagement portion 710 to intersect the proximal surface 712.

When the chamfer-cutting guide 700 is attached to the broach 674, a central axis of the aperture 732 may align with the axis of the femoral neck. The chamfer-cutting guide 700 may engage the broach 674 via a spring loaded mechanism, which may be activated by pressing on a button 765 or lever. When button 765 becomes depressed, the larger aperture 732 of chamfer guide 700 may become able to slide over the proximal portion of the broach. When the external force on the button 765 is released, the spring force may push a side portion of the larger aperture 732 rigidly against the broach post, thus compressing it against the side of the larger

aperture 732. This may rigidly lock the cutting guide 700 into place. The engagement portion 710 may include additional apertures or connecting features to secure the chamfer-cutting guide 700 during use. Examples of additional connecting features may include a protruded portion, a morse-taper or other male, or female, connecting features.

The second step 722 and third step 724 may be connected by a second vertical surface 734. The third step 726 may intersect the second vertical surface 734 and a third vertical surface 736 that is opposite the second vertical surface 734. The third vertical surface 736 may connect to a second rounded corner 738, where it intersects a distal surface 752 of the guide portion 750.

A bottom perspective view of chamfer cutting guide 700 is shown in FIG. 25. The guide portion 750 may be substantially elongated and generally rectangular in shape. The guide portion 750 may include a distal surface 752 that is substantially rectangular. The distal surface 752 may extend between the engagement portion 710 and a flat end portion 715. A cross-section of guide portion 750 is illustrated in FIG. 26. The cross section of guide portion 750 is shown to be generally triangular and contains a central body 754 that may be an equilateral triangle, and two appendages 756, 758 that sit superior and have interior edges 766, 768 that are parallel to the edges of the triangular body 754. The space between the appendages 756, 758 and the edges of the triangular body 754 make up channels or slots 760, 762 that act as cutting guides to establish a predetermined cutting trajectory for a traditional surgical saw or other surgical cutting instrument. The two cutting guides 760, 762 may act as an anterior slot and a posterior slot for creating the desired resulting bone resection on both anterior and posterior sides of the proximal femur. The two slots 760, 762 may intersect at a proximal primary channel 764 that extends along the length (of the zaxis) of the central body 754, as illustrated in FIG. 27A.

The appendages 756, 758 may have exterior surfaces 770, 772 that serve as secondary cutting guides that allow for greater surgeon flexibility. The flat exterior surfaces 770, 772 create additional guiding interface that a powered saw can rest on while creating an alternate bone resection.

Referring to FIGS. 27A-27B a proximal femur is shown with an attached chamfer cutting guide 700 after anterior and posterior resections that were performed using the cutting guide 700. The resultant geometry post-resection creates flat surfaces 787, 789 on anterior/posterior surfaces of the femur. This is similar in nature to the design of a traditional hip stem that includes flat surfaces on the anterior/posterior sides of the neck to avoid impingement with the cup and liner.

Referring to FIG. 28, a diagram illustrates methods of preparing a femur with anterior and/or posterior chamfer cuts. The method starts with step 800, although it is appreciated that this method may occur after another method, such as a method for establishing a surgical approach to a hip joint. In step 810, an initial neck resection is made by removing the femoral head at or near the head-neck junction, or in a subcapital region of the femoral neck. The intramedullary canal is then opened in step 820. In step 830, the proximal metaphysis is then broached. Progressively larger broaches may be used sequentially in step 830. The broach may be positioned in the proximal femur such that it is flush with the resection level. Referring to step 840, one or more chamfer cuts are made at one or more locations (anterior, posterior, medial) around the neck resection. At this point, the method ends with step 850, although it is appreciated that this method may occur before another method, such as a method for implanting a prosthetic hip stem.

In optional step 822, a final neck resection may be made with reference to the final broach in its final position in the proximal femur. Step 822 may occur after step 830.

Referring to optional step 832, the distal intramedullary canal, or diaphysis, may be reamed before or after steps 822 and/or 830. Progressively larger reamers may be used sequentially in step 832.

Steps 842, 844, and 846 describe various possible techniques and instruments that may be used for creating planar chamfer cuts. Referring to step 842, a planar chamfer cut may be made with reference to a guide surface of the broach, similar to the system described in FIGS.21A-21B. Referring to step 844, the planar chamfer cut may alternatively be made with reference to a guide surface of a cutting guide, such as the cutting guide described in FIGS. 24-27B. Referring to step 846, the

planar chamfer cut may also be made with a calcar reamer, such as the calcar reamer described in FIGS. 17A-17C, or with a saw. Referring to step 848, a curved chamfer cut may be made with reference to the broach using a curved chamfer cutting guide surface, such as that described on the instrument illustrated in FIGS. 18A-19B.

The high impact femoral hip stem incorporates neck-strengthening adaptations, such as the shorter neck length and wider neck diameter, which improve the stem's ability to sustain dynamic loads in service. The stem also incorporates stress-transfer adaptations, such as a higher, medialized resection level and porous coating level, which improve the stem's ability to transfer stress to the calcar region of the femoral neck. The stem features a pronounced medial flare which also aids in stress transfer by increasing the projected area of the stem when viewed along the stem axis. Because the stem transfers more stress to the calcar region, the stress distribution in the remainder of the femur may better approximate that of the natural femur, thus preserving more bone substance over a longer period of time.

The reference to any prior art in this specification is not, and should not be taken as, an acknowledgement or any form of suggestion that the prior art forms part of the common general knowledge in Australia.

In this specification, the terms "comprise", "comprises", "comprising" or similar terms are intended to mean a non-exclusive inclusion, such that a system, method or apparatus that comprises a list of elements does not include those elements solely, but may well include other elements not listed.

#### **CLAIMS**

1. A system for creating a bone resection on a proximal femur with a surgical cutting instrument for implanting a femoral hip prosthesis, comprising:

a broach, the broach is shaped to be inserted into a proximal femur to form a hole, the broach including a cutting guide connection feature defined on a surface of the broach; and

a cutting guide configured to be connected to the broach, the cutting guide including an engagement portion and a guide portion, the engagement portion including an underside surface defining a broach connection feature sized and configured to receive and couple to the cutting guide connection feature of the broach the guide portion having an elongated structure that extends from the engagement portion, the elongated structure exhibiting a first elongated slot extending laterally through the elongated structure and along a length of the elongated structure, the elongated structure having a first guide surface extending along the first elongated slot in the elongated structure such that the first guide surface defines a first resection plane acutely angled relative to the underside surface of the cutting guide; wherein the elongated structure longitudinally extends transverse relative to an axis defined by the broach connection feature, wherein the first elongated slot acts as a cutting guide slot configured to receive the surgical cutting instrument for creating the bone resection on the proximal femur.

- 2. The system of claim 1, wherein the elongated structure of the guide portion comprises a central elongated body having a triangular cross-section, the central elongated body defining the first guide surface and a second guide surface..
- 3. The system of claim 2, wherein the guide portion comprises a first elongated appendage extending along the first guide surface to define the first elongated slot therebetween and a second elongated appendage extending along the second guide surface to define a second elongated slot therebetween.

- 4. The system of claim 1, wherein the cutting guide connection feature comprises a post, wherein the broach connection feature comprises a socket, wherein the socket is shaped to receive a portion of the post.
- 5. The system of claim 1, wherein the underside surface of the engagement portion comprises a first step, wherein the first step substantially prevents axial rotation of the cutting guide relative to the broach.
- 6. The system of claim 4, wherein the engagement portion comprises an actuatable button, wherein when an external force is applied to the button, the socket is able to receive a portion of the post, wherein when the external force is released, the engagement portion engages the post such that the post is secured within the socket.
- 7. The system of claim 1, further comprising a prosthesis, wherein a portion of the prosthesis is shaped to fit into the hole.
- 8. A system for implanting a femoral hip prosthesis, comprising:
  - a broach, the broach shaped to be inserted into a proximal femur to form a hole, the broach including a cutting instrument connection feature defined on a surface of the broach; and
  - a cutting instrument configured to be connected to the broach, the cutting instrument including an engagement portion and a bone cutting surface, the bone cutting surface having a cup shaped inner surface extending from the engagement portion, the engagement portion including a broach connection feature at an apex of the cup shaped inner surface, the broach connection feature sized and configured to directly mate to the cutting instrument connection feature of the broach,

wherein the cutting instrument is rotatable such that, upon the instrument being rotated, the cup shaped inner surface of the bone cutting surface engages and removes a section of bone.

- 9. The system of claim 8, wherein the cutting instrument connection feature comprises a proximal socket, wherein the broach connection feature comprises a post, wherein the post is shaped to be received within the proximal socket.
- 10. The system of claim 8, wherein the bone cutting surface is conical.
- 11. The system of claim 8, wherein the cutting instrument comprises a handle, wherein the handle is connected to the engagement portion.
- 12. The system of claim 8, wherein the bone cutting surface comprises at least one aperture, wherein the at least one aperture comprises a cutting edge.
- 13. The system of claim 3, wherein the second guide surface defines a second resection plane acutely angled relative to the underside surface of the cutting guide and the first guide surface.
- 14. The system of claim 1, wherein the elongated structure longitudinally extends transverse relative to an axis defined by the broach connection feature.
- 15. The system of claim 8, wherein the cup shaped inner surface of the bone cutting surface comprises a plurality of apertures each at least partially defined with a cutting edge.
- 16. The system of claim 8, wherein the cutting instrument connection feature comprises a post, wherein the broach connection feature comprises a socket, wherein the post is shaped to be received within the proximal socket.
- 17. The system of claim 8, wherein the cup shaped inner surface extends radially from the apex with varying degrees of curvature.




Fig. 1 (Background Art)




Fig. 2 (Prior Art)




Fig. 3 (Prior Art)




Fig. 4 (Prior Art)




Fig. 5A (Prior Art)

Fig. 5B (Prior Art)




Fig. 6A




Fig. 6B




Fig. 7




Fig. 8




Fig. 9





Fig. 10A



IM Diameter (mm)	Standard Offset (SO) (mm)	High Offset (HO) (mm)	High Impact High Offset (HIHO) (mm)	High Offset Plus (HO <sup>+</sup> ) (HO+5mm)
9	34	40	40	45
10	35	41	41	46
11	36	42	42	47
12	37	43	43	48
13	37	43	43	48
14	38	46	46	51
15	39	47	47	52
16	40	48	48	53
17	40	48	48	53
18	41	49	49	54

Fig. 11B




Fig. 12A




Fig. 12B




Fig. 13A




Fig. 13B




Fig. 13C




Fig. 13D




Fig. 13E






Fig. 14A






Fig. 14B




FIG. 15C




FIG. 15D




FIG. 15A




FIG. 15B





FIG. 16B

FIG. 16A






FIG. 16C




FIG. 16D




FIG. 16E




FIG. 16F

19/30




FIG. 17A




FIG. 17B




FIG. 17C




Fig. 18B




Fig. 19A




Fig. 19B




Fig. 20A




Fig. 20B




FIG. 21A







FIG. 21B










Fig. 24




Fig. 25





Fig. 26






FIG. 28