



(11)

EP 3 033 418 B1

(12)

EUROPEAN PATENT SPECIFICATION

(45) Date of publication and mention
of the grant of the patent:
11.12.2019 Bulletin 2019/50

(51) Int Cl.:
C12N 5/0797^(2010.01) A61K 35/30^(2015.01)

(21) Application number: **14755411.7**

(86) International application number:
PCT/GB2014/052509

(22) Date of filing: **14.08.2014**

(87) International publication number:
WO 2015/022545 (19.02.2015 Gazette 2015/07)

(54) **STEM CELL MICROPARTICLES AND MIRNA**

STAMMZELLENMIKROPARTIKEL UND MIRNA

MICROPARTICULES DE CELLULES SOUCHES ET MIARN

(84) Designated Contracting States:
**AL AT BE BG CH CY CZ DE DK EE ES FI FR GB
GR HR HU IE IS IT LI LT LU LV MC MK MT NL NO
PL PT RO RS SE SI SK SM TR**

(30) Priority: **14.08.2013 GB 201314573**
09.10.2013 GB 201317887

(43) Date of publication of application:
22.06.2016 Bulletin 2016/25

(73) Proprietor: **Reneuron Limited**
Guildford
Surrey GU2 7AF (GB)

(72) Inventors:
• **HICKS, Caroline**
Guildford
Surrey GU2 7AF (GB)
• **SINDEN, John**
Guildford
Surrey GU2 7AF (GB)
• **STEVANATO, Lara**
Guildford
Surrey GU2 7AF (GB)
• **CORTELING, Randolph**
Guildford
Surrey GU2 7AF (GB)

(74) Representative: **Elkington and Fife LLP**
Prospect House
8 Pembroke Road
Sevenoaks, Kent TN13 1XR (GB)

(56) References cited:
WO-A1-2013/150303 WO-A1-2014/013258

WO-A2-03/046141 KR-A- 20100 081 003

- **DUKJIN KANG ET AL: "Proteomic Analysis of Exosomes from Human Neural Stem Cells by Flow Field-Flow Fractionation and Nanoflow Liquid Chromatography-Tandem Mass Spectrometry", JOURNAL OF PROTEOME RESEARCH, vol. 7, no. 8, 1 August 2008 (2008-08-01), pages 3475-3480, XP055066064, ISSN: 1535-3893, DOI: 10.1021/pr800225z**
- **CHENJIE YANG ET AL: "Immunosuppressive Exosomes: A New Approach for Treating Arthritis", INTERNATIONAL JOURNAL OF RHEUMATOLOGY, vol. 176, no. 12, 1 January 2012 (2012-01-01), pages 7385-8, XP055155835, ISSN: 1687-9260, DOI: 10.1038/labinvest.2010.152**
- **E. VAN DER POL ET AL: "Classification, Functions, and Clinical Relevance of Extracellular Vesicles", PHARMACOLOGICAL REVIEWS, vol. 64, no. 3, 1 July 2012 (2012-07-01), pages 676-705, XP055155968, ISSN: 0031-6997, DOI: 10.1124/pr.112.005983**
- **INDIRA VISHNUBHATLA ET AL: "The Development of Stem Cell-derived Exosomes as a Cell-free Regenerative Medicine", JOURNAL OF CIRCULATING BIOMARKERS, vol. 3, no. 2, 30 April 2014 (2014-04-30), pages 1-14, XP055155696, DOI: 10.5772/58597**
- **Valentina Fonsato ET AL: "Human Liver Stem Cell-Derived Microvesicles Inhibit Hepatoma Growth in SCID Mice by Delivering Antitumor MicroRNAs", Stem Cells, vol. 30, no. 9, 20 September 2012 (2012-09-20), pages 1985-1998, XP055161662, ISSN: 1066-5099, DOI: 10.1002/stem.1161**

Note: Within nine months of the publication of the mention of the grant of the European patent in the European Patent Bulletin, any person may give notice to the European Patent Office of opposition to that patent, in accordance with the Implementing Regulations. Notice of opposition shall not be deemed to have been filed until the opposition fee has been paid. (Art. 99(1) European Patent Convention).

EP 3 033 418 B1

- Stefania Bruno ET AL: "Effects of Mesenchymal Stromal Cell-Derived Extracellular Vesicles on Tumor Growth", *Frontiers in Immunology*, vol. 5, 11 August 2014 (2014-08-11), XP055161795, DOI: 10.3389/fimmu.2014.00382

Description**Field of the Invention**

5 **[0001]** This invention relates to neural stem cell microparticles, their use and production thereof, in particular neural stem cell microparticles and their use in therapy of glioma, in particular glioblastoma.

Background of the Invention

10 **[0002]** Stem cells have the ability to self-renew and to differentiate into functionally different cell types. They have the potential to be a powerful therapeutic tool, for example in the growing field of Regenerative Medicine, in particular regenerative therapy requiring tissue replacement, regeneration or repair (Banerjee et al. 2011). Endogenous stem cells have also been implicated as targets (endogenous "cancer stem cells") of anti-cancer therapy, where it is proposed to treat the cancer by eliminating the cancer stem cells that are thought to drive cancer growth and metastasis. More recently, engineered mesenchymal stem cells have been proposed as delivery vehicles in anti-cancer therapy (Dai et al., 2011; Shah et al. 2012). However, there are drawbacks to the use of stem cells in therapy: there is a need for a consistent and substantial supply of stem cells with functional and phenotypic stability and the associated high costs and time delay caused by cell generation, storage, transport and handling; there is a requirement for immunological compatibility to avoid rejection of the stem cells by the recipient; and there are complex regulatory issues related to potential safety risks of tumour or ectopic tissue formation. Further, despite the therapeutic efficacy of stem cell transplantation, there is no convincing evidence for a direct long-term effect of the transplanted stem cells, for example through engraftment and differentiation into reparative or replacement cells.

25 **[0003]** Neural stem cells (NSCs) are self-renewing, multipotent stem cells that generate neurons, astrocytes and oligodendrocytes (Kornblum, 2007). The medical potential of neural stem cells is well-documented. Damaged central nervous system (CNS) tissue has very limited regenerative capacity so that loss of neurological function is often chronic and progressive. Neural stem cells (NSCs) have shown promising results in stem cell-based therapy of neurological injury or disease (Einstein et al. 2008). Implanting neural stem cells (NSCs) into the brains of post-stroke animals has been shown to be followed by significant recovery in motor and cognitive tests (Stroemer et al. 2009). It is not completely understood how NSCs are able to restore function in damaged tissues but it is now becoming increasingly recognised that NSCs have multimodal repairing properties, including site-appropriate cell differentiation, pro-angiogenic and neurotrophic activity and immunomodulation promoting tissue repair by the native immune system and other host cells (Miljan & Sinden, 2009, Horie et al., 2011). It is likely that many of these effects are dependent on transient signalling from implanted neural stem cells to the host milieu, for example NSCs transiently express proinflammatory markers when implanted in ischaemic muscle tissue damage which directs and amplifies the natural pro-angiogenic and regulatory immune response to promote healing and repair (Katare et al., Clinical-grade human neural stem cells promote reparative neovascularization in mouse models of hindlimb ischemia. Arteriosclerosis, Thrombosis and Vascular Biology, in press). In chronic stroke brain, NSCs also have a substantial neurotrophic effect. For example, they promote the repopulation of the stroke-damaged striatal brain tissue with host brain derived doublecortin positive neuroblasts (Hassani, O'Reilly, Pearse, Stroemer et al., PLoS One. 2012;7(11)).

40 **[0004]** Furthermore, on the basis of a large body of NSC restorative effects in animal models with chronic stroke, a clinical trial using neural stem cells is being carried out by ReNeuron Limited (Surrey, UK), to trial the treatment of disabled stroke patients using its "CTX0E03" conditionally-immortalised cortex-derived neural stem cells (Clinicaltrials.gov Identifier: NCT01151124).

45 **[0005]** Mesenchymal stem cells (MSCs) are lineage-restricted stem cells which have the potential to differentiate into mesenchymal cell types only, namely of the adipocytic, chondrocytic and osteocytic lineages (Pittenger et al. 1999; Ding et al. 2011). MSCs (also referred to as Mesenchymal Stromal Cells and Mesenchymal Progenitor Cells) are derived from a variety of sources including bone marrow, blood, adipose and other somatic tissues. The therapeutic potential of MSCs, however, is more directed towards the application of their pro-angiogenic and immune modulating properties as undifferentiated cells. Production of human MSCs is limited by the inability of these cells to expand in numbers stably beyond approximately 15-20 population doublings.

50 **[0006]** Mesenchymal stem cell-conditioned medium (MSC-CM) has a therapeutic efficacy similar to that of MSCs themselves, suggesting a paracrine mechanism of MSC-based therapy (Timmers et al. 2007). WO-A-2009/105044 discloses that particles known as exosomes, secreted by MSCs, comprise at least one biological property of the MSCs and suggests the use of these MSC particles in therapy, while Théry et al. 2011 provides a general review of exosomes and other similar secreted vesicles. Whereas some of the drawbacks of using stem cells directly as therapeutic agents are overcome by using the mesenchymal stem cell-derived exosomes (e.g. storage, transport and handling), the problem remains of providing a consistent and substantial supply of functionally and phenotypically stable stem cells to produce the exosomes. For therapeutic use, the exosomes preferably need to be produced on a large scale. In the absence of

a stem cell line, replenishment of the cells through repeated derivation from a source of stem cells is required, which incurs recurring costs for testing and validation of each new batch. Furthermore, the diseases and disorders that can be treated by MSCs may be limited.

[0007] WO-A-2013/150303 and WO-A-2014/013258 disclose microparticles produced by neural stem cells, methods for making those microparticles and uses of those microparticles, in particular for use in regenerative therapy. Fonsato et al (Stem Cells. 2012 Sep;30(9):1985-98) report that microvesicles derived from human adult liver stem cells (HLSC) may reprogram *in vitro* HepG2 hepatoma and primary hepatocellular carcinoma cells by inhibiting their growth and survival.

[0008] There remains a need for improved stem cell-based therapies.

Summary of the Invention

[0009] The invention provides a neural stem cell exosome for use in a method of treating glioma, in particular glioblastoma.

[0010] The disclosure is based on the surprising finding that neural stem cells contain microparticles that are therapeutically useful. In particular, the inventors have surprisingly identified neural stem cell microparticles that are able to: inhibit cell migration of fibroblasts; inhibit migration of cancer cells; induce differentiation of cancer cells; and/or induce or enhance an immune response against cancer cells. These properties make the neural stem cell microparticles suitable for use in therapy, in particular for treating cancer.

[0011] Cell migration is well-known to play an important role in the progression of diseases such as cancer (for example during angiogenesis, tumour formation, metastasis and tissue invasion), fibrosis (for example during the accumulation of fibroblasts in the fibrotic tissue), atherosclerosis and rheumatoid arthritis. Microparticles that inhibit cell migration are therefore useful in the treatment or prevention of diseases that involve unwanted cell migration, such as cancer, in particular metastatic cancer, fibrosis, atherosclerosis and rheumatoid arthritis.

[0012] Microparticles of the disclosure are shown, in the Examples, to inhibit fibroblast migration. Fibroblasts and the migration of fibroblasts are known to play a role in angiogenesis and so the microparticles of the disclosure, which inhibit fibroblast migration, are also useful for use in the therapy of unwanted or undesirable angiogenesis.

[0013] Additionally, the Examples show that glioblastoma cells, pre-treated *in vitro* for 24 hours with neural stem cell exosomes did not engraft into the striatum of Balb-C mice *in vivo*. Histopathology demonstrated the presence of necrotic cell bodies at the site of implantation and evidence of a host cellular response. These data indicate that these microparticles are suitable for use in the treatment of cancer, particularly a cancer of the CNS such as a glioblastoma, by promoting the destruction of cancer cells by the immune system.

[0014] Neural stem cell microparticles that are able to inhibit fibroblast cell migration and induce or enhance an immune response against cancer cells are also disclosed herein and have been isolated from neural stem cells cultured in a multi-compartment bioreactor for 11 weeks. Accordingly, one way to obtain these neural stem cell microparticles is to isolate them from neural stem cells that have been cultured in a multi-compartment bioreactor for at least 10 weeks, for example 71 days or more. The microparticles of the disclosure may also be obtained from other culture conditions and periods, in particular culture conditions that allow stem cell differentiation.

[0015] The Examples further show that tumour (glioblastoma U373) cells show significantly reduced migration when treated with neural stem cell exosomes. The exosomes of the invention may therefore be used to treat a glioblastoma, by inhibiting tumour cell migration.

[0016] Additionally, neural stem cell exosomes are shown in the Examples to promote differentiation of tumour (glioblastoma U373) cells *in vitro*. The Examples also show this differentiation *in vivo*, where tumour (glioblastoma U373) cells, treated with neural stem cell exosomes and implanted into mouse brains, demonstrate a reduction in the stem cell marker nestin. Cancer stem cells drive tumourigenesis, are linked with metastasis, high grade and poor prognosis. A more differentiated tumour typically correlates with improved prognosis, so exosomes that are able to effect differentiation are expected to be useful in the treatment of cancer. Therefore, the ability of microparticles isolated from neural stem cells to reduce the stemness of cancer cells indicates that these microparticles are useful in the treatment of cancer, in particular a cancer that is positive for nestin expression such as melanoma, breast cancer or glioblastoma. Typically, the cancer is a cancer of the CNS such as a glioblastoma. Nestin is reported to correlate with aggressive growth, metastasis, and poor prognosis in cancers, so agents that reduce nestin expression are greatly needed. Neural stem cell microparticles that are able to inhibit tumour cell migration and promote differentiation of tumour cells have been isolated from a neural stem cell line cultured under standard conditions. Accordingly, one way to obtain neural stem cell microparticles that are able to inhibit tumour cell migration and promote differentiation of tumour cells is to isolate them from neural stem cells that have been cultured under standard conditions. These cells may be from the CTX0E03 cell line (deposited with the ECACC as Accession No. 04091601). The standard culture conditions typically maintain the characteristics of the cell line, in particular the stemness of the cell line, typically do not permit differentiation, and typically provides proliferating cells. Typically, the cells proliferate with a doubling time of 2 to 4 days and are passaged when

sub-confluent.

[0017] The Examples include a pilot *in vivo* study of the administration of exosomes of the invention to human glioblastoma xenografts, observing tumour sensitivity to the microparticles, a trend towards a reduction in tumour volume, and increased survival. Histopathology of the tumour cells shows, in one animal, a particularly dramatic and effective ablation of the tumour mass. The Examples also provide Next Generation Sequence (NGS) analysis of the miRNA content of neural stem cell exosomes. One of the Examples revealed the presence of a set of miRNAs: hsa-mir-1246, hsa-mir-4488, hsa-mir-4492, and hsa-mir-4532, each of which is shown to reduce glioma cell proliferation. These data provide further evidence that microparticles containing the miRNAs may be used to treat cancer.

[0018] A first aspect of the invention provides a neural stem cell exosome for use in a method of treating glioma. The glioma may be glioblastoma. The neural stem cell exosome may inhibit glioblastoma cell migration; and/or induce differentiation of a glioblastoma cell, typically a glioblastoma cell that is positive for nestin expression. In another aspect described herein, the exosome promotes destruction of glioblastoma tumour cells by inducing or enhancing an immune response against the tumour cells.

[0019] A microparticle may be an exosome, microvesicle, membrane particle, membrane vesicle, exosome-like vesicle, ectosome-like vesicle, ectosome or exovesicle. According to the invention, the microparticle is an exosome. The microparticle disclosed herein may be derived from a neural stem cell that has been cultured in an environment that allows stem cell differentiation. The microparticle disclosed herein may or may not be isolated from partially-differentiated neural stem cells; as discussed below, the presence of GFAP (an astrocyte marker) or DCX (an early neuronal marker) on the cells indicates that the neural stem cells have begun to differentiate. As disclosed herein, an environment that allows stem cell differentiation is a multi-compartment bioreactor. The microparticle disclosed herein may be isolated from neural stem cells that have been cultured in a multi-compartment bioreactor for at least 10 weeks. The microparticle may be isolated from cultured neural stem cells that have been confluent on the membrane of a multi-compartment bioreactor for at least one week, at least 2 weeks, typically at least 3 weeks, at least 4 weeks, at least 5 weeks or more. Conversely, microparticles of the invention can be produced from neural stem cells that have not begun to differentiate, for example by isolation from sub-confluent cultured neural stem cells, or by isolation from cells that have been confluent for less than one week on the membrane of a multi-compartment bioreactor or in a standard cell culture flask such as a T-175 flask. As used herein, the term "confluent" is given its usual meaning in the art, wherein the cells in the culture are all in contact and have no further room to grow; confluent cells cover substantially all of the membrane in the multi-compartment bioreactor.

[0020] The microparticle may be derived from a neural stem cell line. In some embodiments, the neural stem cell line may be the "CTX0E03" cell line, the "STR0C05" cell line, the "HPC0A07" cell line or the neural stem cell line disclosed in Miljan et al Stem Cells Dev. 2009. In some embodiments, the microparticle is derived from a stem cell line that does not require serum to be maintained in culture. The microparticle may have a size of between 30 nm and 1000 nm, or between 30 and 200 nm, or between 30 and 100 nm, as determined by electron microscopy; and/or a density in sucrose of 1.1-1.2 g/ml. The microparticle may comprise RNA. The RNA may be mRNA, miRNA, and/or any other small RNA. The microparticle may comprise one, two, three or four of hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532; alternatively, it may comprise 1, 2, 3, 4 or 5 of hsa-miR-181a-5p, hsa-miR-1246, hsa-miR-127-3p, hsa-miR-21-5p and hsa-miR-100-5p; or it may comprise 1, 2, 3, 4 or 5 of hsa-miR-181a-5p, hsa-let-7a-5p, hsa-let-7f-5p, hsa-miR-92b-3p, and hsa-miR-9-5p. The microparticle may comprise one or more lipids, typically selected from ceramide, cholesterol, sphingomyelin, phosphatidylserine, phosphatidylinositol, phosphatidylcholine. The microparticle may comprise one or more tetraspanins, typically CD63, CD81, CD9, CD53, CD82 and/or CD37. The microparticle may comprise one or more of TSG101, Alix, CD109, thy-1 and CD133. The microparticle may comprise at least 10 of the proteins present in Table 20 or Table 22. The microparticle may comprise at least one biological activity of a neural stem cell or a neural stem cell-conditioned medium. At least one biological activity may be an anti-cell migration activity, a pro-differentiation activity or an anti-angiogenic activity. The microparticle of the invention is typically isolated or purified.

[0021] The disease to be treated is a cancer of the CNS, namely a glioma. An exemplary glioma is a glioblastoma, which may be a giant cell glioblastoma or a gliosarcoma.

[0022] The neural stem cell exosome is used to treat the cancer.

[0023] In a further embodiment, the exosomes of the invention treat the cancer by inhibiting migration of the cancer cells.

[0024] In yet a further embodiment, the exosomes of the invention treat the cancer by inducing differentiation of the cancer cells. Typically, differentiation is induced in cancer cells that express nestin.

[0025] As described, the exosomes disclosed herein treat the cancer by inducing or enhancing an immune response against the cancer cells. The immune response typically comprises the activation and/or proliferation of glial cells such as microglia.

[0026] In one embodiment, the therapeutic exosome is isolated from proliferating neural stem cells that have been cultured under conditions that typically maintain the characteristics of the cell line, in particular the stemness of the cell line. These are typically the standard culture conditions for a given cell or cell line, which do not permit differentiation of the stem cells. Typically, proliferating cells have a doubling time of 2 to 4 days. These neural stem cells are typically

passed when sub-confluent.

[0027] As disclosed herein, it has also been found that it is possible to alter the production of microparticles by stem cells, by culturing the stem cells (optionally for at least 10 weeks) and adding components to the culture medium, by culturing the stem cells (optionally for at least 10 weeks) under hypoxic conditions, or by co-culture with other cell types (optionally for at least ten weeks).

[0028] Accordingly, the disclosure provides a method of producing a stem cell microparticle that inhibits cell migration, typically a neural stem cell microparticle that inhibits cell migration. The stem cells may be cultured under conditions that allow the efficient removal of metabolic waste. The method may comprise culturing the stem cells for at least 10 weeks in an environment that allows stem cell differentiation and collecting the microparticles that are produced by the cells; microparticles produced by this method are typically able to inhibit fibroblast migration. The microparticles may be isolated from partially-differentiated neural stem cells. In one embodiment, an environment that allows stem cell differentiation is culture in a multi-compartment bioreactor, typically for a prolonged period of time, for example more than seven days and usually more than ten weeks.

[0029] The method may alternatively comprise culturing the cells under conditions that do not allow differentiation to occur, and collecting the microparticles that are produced by the cells; microparticles produced by this method are typically able to inhibit glioblastoma migration.

[0030] The method may comprise isolating a microparticle from a stem cell-conditioned medium. The stem cell-conditioned medium may comprise one or more additive components or agents which stimulate the release of microparticles by the stem cells into the medium. The one or more components may be selected from transforming growth factor-beta (TGF- β), interferon-gamma (IFN- γ) and/or tumour necrosis factor-alpha (TNF- α). The microparticles may be isolated from stem cell-conditioned medium wherein the stem cells were cultured under hypoxic conditions. The microparticles may be isolated from stem cell-conditioned medium produced by stem cells co-cultured with a different cell type, typically endothelial cells, in order to create the NSC niche environment.

[0031] The disclosure also provides a composition comprising a neural stem cell microparticle according to the first aspect and a pharmaceutically acceptable excipient, carrier or diluent. In one embodiment, the microparticle of the invention inhibits glioblastoma cell migration, typically as determined in a transmembrane or wound healing (scratch) assay. In another embodiment, the microparticle of the invention induces differentiation of a glioblastoma cell, typically as determined by cell morphology and/or marker expression. A decrease in the stem cell marker nestin typically indicates differentiation.

Brief Description of the Drawings

[0032]

Figure 1 shows the effect of neural stem cell exosome treatment on human dermal fibroblast migration in a transmembrane assay. The top panel depicts the assay apparatus and the bottom panel compares the number of cells that migrated through the membrane in the presence of medium alone (a, "basal"), in the presence of 20 μ g/ml exosomes isolated from "0"-week CTX0E03 neural stem cells (b, "Exosome (0)") and in the presence of 20 μ g/ml exosomes isolated from CTX0E03 neural stem cells cultured for 11 weeks in the Integra CELLLine AD1000 culture system (c, "Exosome (11)"), determined after 6 hours and after 24 hours assay incubation.

Figure 2 depicts the human dermal fibroblast cells that migrated through the membrane in the presence of each of the basal, 0-week and 11-week exosomes.

Figure 3A shows the results of a wound closure/scratch assay representing the migration activity of normal human dermal fibroblasts (NHDF) in response to conditioned medium from CTX0E03 cells cultured for 2 weeks and exosomes purified from the conditioned medium of CTX0E03 cultured for 2 weeks in the Integra CELLLine AD1000 culture system. Figure 3B shows the results of a scratch assay after 72 hours, comparing the effect of 10 μ g 2-week CTX0E03 exosomes to basal conditions (without exosomes). Figure 3C shows the % of healed areas for basal conditions, 2 μ g/ml exosomes, 6 μ g/ml exosomes, 20 μ g/ml exosomes and an LSGS (low serum growth supplement) positive control. The top panel of Figure 3C shows exosomes isolated from CTX0E03 cells cultured for 2 weeks in the Integra Celline system and the bottom panel of Figure 3C shows exosomes isolated from CTX0E03 cells cultured for 6 weeks in the Integra Celline system. Figure 3D compares 2-week CTX0E03 cells to a negative control (saline) in an *in vivo* injection wound healing assay.

Figure 4 depicts electron micrographs of CTX0E03 conditionally-immortalised neural stem cells producing microparticles. Panels A-E show intracellular multivesicular bodies (MVBs) containing exosomes between 30nm and 50nm in diameter and Panel F shows microvesicles >100nm in diameter released from neural stem cells through

a process of budding at the cell membrane.

Figure 5 is an outline protocol for the identification, characterisation and production of microparticles from stem cells.

Figure 6 shows the FACS detection (at 2ug/ml, 1:250) of (i) CD63 in 2-week Integra cultured CTX0E03 exosomes (top left panel) and microvesicles (top right panel) and (ii) CD81 in 2-week Integra cultured CTX0E03 exosomes (bottom left panel) and microvesicles (bottom right panel).

Figure 7 shows the results of NanoSight analysis undertaken to determine the particle size and concentration of CTX0E03 exosomes (Figure 7A) and microvesicles (Figure 7B) cultured in the Integra Celline system for 1, 2, 3, 4, 5 and 6 weeks.

Figure 8A shows the amount of protein (measured by BCA assay) extracted from 15ml of media containing microparticles purified from the Integra system compared to normal culture conditions (3 days T175). Figure 8B shows the amount of isolated total RNA measured at 260/280nm extracted from 15ml of CTX0E03 conditioned media containing microparticles purified by filtration from the Integra system compared to normal culture conditions (3 days T175).

Figure 9 shows the quantity of purified exosomes obtained per ml culture medium from standard CTX0E03 (T175) cultures vs. the Integra CELLline system at the 3 week time point.

Figure 10A shows the concentration of exosomes harvested from two different flasks after 1 week, 2 weeks and 3 weeks of CTX0E03 Integra CELLline culture system. Figure 10B shows the concentration of exosomes harvested from a single Integra CELLline flask during a 6 week continuous culture of CTX0E03 cells.

Figure 11 shows the fold change of expression levels of various mRNA markers measured in CTX0E03 cells cultured for 3 weeks in the Integra CELLline system compared to standard ("control") CTX0E03 (T175) cultures.

Figure 12 shows the fold up and down regulation of various miRNAs in exosomes obtained from CTX0E03 cells cultured for 3 weeks in Integra bioreactor culture and microparticles obtained from standard CTX0E03 (T175) cultures, assessed against a baseline expression level in CTX0E03 cells in standard (T175) culture.

Figure 13 depicts miRNA deep sequencing results. The miRNA profiles obtained from deep sequencing of miRNA from CTX0E03 cells ("CTX"), microvesicles ("MV") and exosomes ("EXO") cultured under standard (T175) conditions are shown in Figure 13A and 13B (results from two standard cultures, "EH" and "EL"). Figure 13C shows the percentage of miRNAs that are up-shuttled, the same, or down-shuttled in the exosomes compared to producer cells, for (i) the standard culture, (ii) 6 week Integra bioreactor culture and (iii) 11 week bioreactor culture (3 samples). Up-shuttled > 2 , same < 2 , and down-regulated < 2 fold change (log2) accordingly. Figures 13D to 13H show the miRNAs that are shuttled into exosomes compared with the cells producing them. Up-shuttled miRNAs are expressed as fold change calculated using the log2 of the normalized ratio of exosomes/cell producer. The normalization is obtained by dividing reads of each miRNA by total miRNA reads. (D) summarises the most abundant miRNAs in exosomes obtained from the standard CTX0E03 cultures ("EH" and "EL"); (E) shows exosomes obtained from CTX0E03 cells cultured for 6 weeks in an Integra bioreactor, and lists up-shuttled miRNAs with more than 250 reads per exosome sample; (F) shows the miRNAs up-shuttled in exosomes when compared with the producer cells cultured for 11 weeks in an Integra bioreactor. 9 miRNA species are up-shuttled, all of which have more than 250 reads; (G) shows a second sample of the miRNAs up-shuttled in exosomes when compared with the producer cells cultured for 11 weeks in an Integra bioreactor. The diagram lists up-shuttled miRNAs with more than 250 reads per exosome sample; and (H) shows a third sample of the miRNAs up-shuttled in exosomes when compared with cell producer cultured for 11 weeks in an Integra bioreactor, showing up-shuttled miRNAs with more than 250 reads per exosome sample.

Figure 14 is an electropherogram showing the total RNA content profile in 2-week CTX0E03 cells, exosomes and microvesicles as determined by Agilent RNA bioanalyser.

Figure 15 is a schematic presentation of the percentage of coding genes fully overlapping exon, and non-coding transcripts located with intron or intergenic sequences (produced by running NGS BAM files against GENCODE sequence data set).

Figure 16 depicts the top ranking preferentially shuttled novel miRNAs in exosomes and MV compared to CTX0E03 producer cells.

Figure 17 shows Venn diagrams comparing the proteomic data from CTX0E03 exosomes and microvesicles (17A and 17B), and comparing neural stem cell exosomes with mesenchymal stem cell exosomes (17C and 17D). Figure 17A illustrates the number of unique proteins within CTX0E03 exosomes and microvesicles, isolated from week 2 Integra culture system. Figure 17B compares the biological processes associated with the identified proteins within the CTX0E03 exosomes and microvesicles. Figure 17C compares the CTX0E03 neural stem cell exosome proteome to a Mesenchymal Stem Cell exosome, and Figure 17D compares biological processes associated with the identified proteins in the MSC derived exosomes with the neural stem cell derived exosomes.

Figure 18 shows the 30 biological processes found to be associated with NSC derived exosomes and not mesenchymal stem cell exosomes.

Figure 19 shows the presence of necrotic cell bodies and evidence of a host cellular response in the striatum of Balb-C mice 24 hours after implantation of glioblastoma U373 cells that were pre-treated for 24 hours with exosomes isolated from CTX0E03 cells cultured for 11 weeks in a multi-compartment bioreactor.

Figure 20 shows a reduction in nestin expression in glioblastoma U373 cells that have pre-treated in vitro for 24 hours with exosomes isolated from a proliferating culture of CTX0E03 cells and implanted into the striatum of Balb-C mice.

Figure 21 shows that glioblastoma U373 cells that have been treated *in vitro* with exosomes isolated from a proliferating culture of CTX0E03 cells, appear morphologically differentiated and express Glial fibrillary acidic protein (GFAP).

Figure 22 shows that seeding glioblastoma cells together with 20 μ g/ml CTX0E03 exosomes (Figures 22A and 22C) or have been pre-treating glioblastoma cells with 10 μ g/ml CTX0E03 exosomes for 24hours (Figure 22B) reduces glioblastoma migration towards 10% FBS.

Figure 23 shows the inhibitory effects of individual miRNAs on the proliferation of glioma cells: (A) plot of percentage of U373MG cell proliferation, compared to 24hrs control, measured by CyQUANT assay following transfection with hsa-mir-1246, hsa-mir-4488, hsa-mir-4492, or hsa-mir-4532; (B) plot of percentage of U373MG cell proliferation, compared to 0hr control, measured by CyQUANT assay following transfection with hsa-mir-1246, hsa-mir-4488, hsa-mir-4492, or hsa-mir-4532; and (C) plot of percentage of U87 cell proliferation, compared to 0hr control, measured by CyQUANT assay following transfection with hsa-mir-1246, hsa-mir-4488, hsa-mir-4492, or hsa-mir-4532.

Figure 24 shows glioblastoma xenograft individual tumour volumes of mice on the day of assignment to each treatment group.

Figure 25 shows the mean body weights of the mice during the xenograft study. The dotted vertical line indicates the commencement of the dosing phase (on day 12).

Figure 26 summarises the mean tumour volume for the treatment groups measured during the study.

Figure 27 displays the tumour volume data (% pre-dose) of Figure 26 in a truncated format up to study day 25.

Figure 28 shows the final tumour weights, expressed as group mean + standard error of the mean (tumour weight).

Figure 29 shows survival analysis utilising mean tumour diameter (15mm) as the humane survival endpoint.

Figure 30 shows the absolute individual body weights of each mouse in the study.

Figure 31 shows the relative individual body weights of each mouse in the study.

Figure 32 shows the raw data for individual tumour volume measurements.

Figure 33 shows the individual tumour volume plots.

Figure 34 details the tumour weights.

Detailed Description of the Invention

[0033] The invention provides a neural stem cell exosome for use in a method of treating a glioma, for example a glioblastoma.

[0034] The disclosure provides that the present inventors have surprisingly identified that neural stem cells produce microparticles that inhibit cell migration of fibroblasts and cancer cells, induce differentiation of cancer cells, and/or induce or enhance an immune response against the cancer cells. These microparticles are shown to inhibit cell migration and are therefore useful in therapy of diseases comprising unwanted, undesired or deleterious cell migration. The microparticles inhibit fibroblast migration and are therefore also useful in therapy of diseases comprising unwanted, undesired or deleterious angiogenesis, in which fibroblasts play a key role. The microparticles also inhibit tumour cell migration, induce differentiation of tumour cells and enhance an immune response against cancer cells, and are therefore useful in the treatment of cancer. The microparticles of the disclosure can be characterised and identified by these properties, using the assays described herein or other assays known to the skilled person.

[0035] The microparticles are advantageous over the corresponding stem cells because they are smaller and less complex, thereby being easier to produce, maintain, store and transport, and have the potential to avoid some of the regulatory issues that surround stem cells. The microparticles can be produced continuously, by isolation from conditioned media, for example in a bioreactor such as a multi-compartment bioreactor, which allows for large scale production and the provision of an "off-the-shelf" therapy. The multi-compartment bioreactor is typically a two-compartment bioreactor. An exemplary multi-compartment bioreactor is the CeLLine AD1000 bioreactor that is commercially available from Integra Biosciences AG, Zizers, Switzerland (Item No. 90025).

[0036] The inventors have found that the properties of neural stem cell microparticles differ depending on the culture conditions of the stem cells that produce the microparticles, in particular the length of time that the neural stem cells are cultured before the microparticles are harvested. In particular, the inventors have surprisingly identified neural stem cell microparticles that inhibit cell migration and/or induce differentiation of a stem or cancer cell.

[0037] In one disclosed embodiment, microparticles that inhibit fibroblast migration can be isolated from neural stem cells that have been cultured in a multi-compartment bioreactor for at least 10 weeks, e.g. more than 10 weeks. This is particularly surprising because microparticles isolated from the same neural stem cells that have been cultured for less than 10 weeks, for example about 2-6 weeks, have been shown to enhance fibroblast cell migration as seen in wound healing assays.

[0038] In another embodiment, microparticles that are able to induce or enhance a beneficial immune response against cancer cells can also be isolated from neural stem cells that have been cultured in a multi-compartment bioreactor for at least 10 weeks, e.g. more than 10 weeks.

[0039] In a further embodiment, microparticles that are able to reduce tumour cell migration and/or induce cancer or stem cell differentiation are isolated from a proliferating neural stem cell culture. This culture may be in a standard cell culture flask (such as a T-175 flask) or may be in a multi-compartment bioreactor. When the cells producing microparticles of this embodiment are cultured in a multi-compartment bioreactor, they are typically cultured for 4 weeks or less, for example 3 weeks or less, 2 weeks or less, or 1 week or less. This is because, as described elsewhere herein, prolonged culture in a multi-compartment bioreactor allows the stem cells to begin to differentiate, i.e. to express markers for defined neural cell types. Typically, the microparticles that are able to reduce tumour cell migration and/or induce differentiation are isolated from neural stem cells that are negative for markers of differentiated neural cells (e.g. GFAP⁻ and/or DCX⁻) but are positive for one or more markers of neural stem cells (e.g. Nestin⁺).

[0040] Figure 1 (lower panel) and Figure 2 show that exosomes isolated from a non-proliferating CTX0E03 culture significantly abrogate migration of human dermal fibroblasts. This is in contrast to exosomes isolated from a proliferating CTX0E03 culture, which significantly promote migration of human dermal fibroblasts. Accordingly, in one embodiment, microparticles that inhibit cell (e.g. fibroblast) migration may be isolated from non-proliferating neural stem cells. Optionally, these non-proliferating stem cells may be partly differentiated, i.e. express one or more early markers of differentiation. In one embodiment, the neural stem cells from which these microparticles are isolated are positive for DCX (doublecortin), which is an early neuronal marker. In another embodiment, the neural stem cells from which the microparticles are isolated are positive for GFAP (Glial fibrillary acidic protein), which is an astrocyte marker.

[0041] Figure 22 shows that exosomes isolated from a proliferating CTX0E03 culture inhibit the migration of glioblastoma cells towards a positive chemoattractant. Accordingly, in one embodiment, microparticles that inhibit cell (e.g. glioblastoma) migration may be isolated from proliferating neural stem cells that are typically negative for markers of differentiation. In one embodiment, the neural stem cells from which these microparticles are isolated are negative for DCX (doublecortin), which is an early neuronal marker. In another embodiment, the neural stem cells from which the microparticles are isolated are negative for GFAP (Glial fibrillary acidic protein), which is an astrocyte marker.

[0042] Cell migration is well-known to play an important role in the progression of diseases such as cancer (for example

during angiogenesis, tumour formation, metastasis and tissue invasion), fibrosis (for example during the accumulation of fibroblasts in the fibrotic tissue), atherosclerosis and rheumatoid arthritis. The identification of microparticles that are able to inhibit these processes therefore provides a new therapy for these diseases.

[0043] Transmembrane and wound healing assays are physiologically relevant cell-based assays that are predictive of *in vivo* mechanisms of cell migration and allow the identification of compounds that are effective in promoting or inhibiting cell migration, and the recognition of potential undesirable effects. Furthermore there is a good correlation between the results obtained in scratch assays and transmembrane assays (Hulkower et al. 2011), so that these assays can be compared.

[0044] The data presented below demonstrate that microparticles that inhibit cell migration can be isolated from neural stem cells. The microparticles of the disclosure can be produced by any method, not limited to those disclosed or exemplified herein. Whether or not a microparticle is able to inhibit cell migration can be readily determined using the assays described herein.

[0045] It has further been found that, surprisingly, culturing stem cells (of any type, not limited to neural stem cells) in an environment that allows the stem cells to begin to differentiate, increases dramatically the yield of microparticles produced. Typically, the stem cells are NSCs, for example CTX0E03, cultured for at least 10 weeks, for example for 11 weeks, but optionally no more than 20 weeks, 30 weeks or 40 weeks.

[0046] It has also been surprisingly observed that culturing stem cells (of any type, not limited to neural stem cells) in a multi-compartment bioreactor results in partial differentiation of the stem cells, into stem cells in a more differentiated form. This differentiation in culture does not require the addition of an agent to induce differentiation. This differentiation typically requires a culture period of at least one week, at least two weeks, at least three weeks, at least six weeks, at least eight weeks, or at least ten weeks, for example about 11 weeks, but optionally no more than 20 weeks. The changes to the stem cells that occur in culture in a multi-compartment bioreactor are reflected by the microparticles produced by the cultured stem cells. Therefore, by culturing stem cells in a multi-compartment bioreactor, it is possible to induce differentiation of the cells. Accordingly, microparticles from partially differentiated stem cells can be produced by harvesting microparticles from stem cells, for example NSCs such as CTX0E03, cultured in a multi-compartment bioreactor, typically for at least one week, at least two weeks, at least three weeks, at least four weeks, at least five weeks or at least six weeks, at least eight weeks, or at least ten weeks, for example about 11 weeks, but optionally no more than 20 weeks. Typically, the NSCs have been cultured for more than ten weeks. In one embodiment, the disclosure provides a method of producing microparticles by isolating the microparticles from partially-differentiated neural stem cells as described above.

[0047] It has also been found that it is possible to induce the secretion of microparticles from stem cells. Typically, the stem cells are NSCs, for example CTX0E03, typically cultured for at least 10 weeks, for example for 11 weeks or more, but optionally no more than 20 weeks. This finding, which also is not limited to neural stem cells and can be used for the production of microparticles from any stem cell, allows for an improved yield of microparticles to be obtained from a stem cell culture. Several agents have been identified that enhance the secretion of microparticles to different degrees, which has the further advantage of being able to control the amount of microparticles that are secreted. Culturing stem cells under hypoxic conditions also improves microparticle production. Further, it has been found that co-culturing a stem cell with a different cell type, in particular an endothelial cell type can beneficially alter the microparticles that are produced by the stem cell.

[0048] In a further embodiment, the disclosure provides microparticles, typically exosomes, produced by serum-free stem cells. Typically, the stem cells are NSCs, for example CTX0E03, cultured for at least 10 weeks, for example for 11 weeks, but optionally no more than 20 weeks. Serum is required for the successful culture of many cell lines, but contains many contaminants including its own exosomes. As described below, the inventors have produced microparticles from stem cells that do not require serum for successful culture.

Neural Stem Cell Microparticles

[0049] The disclosure provides, in one aspect, microparticles that inhibit cell migration and/or induce differentiation of a stem or cancer cell, obtainable from a neural stem cell. The microparticle is, in one embodiment, obtainable from a neural stem cell that has been cultured in a multi-compartment bioreactor, typically for at least 10 weeks, for example 11 weeks or more, or 12 weeks or more. In another embodiment, the microparticle is obtainable from a proliferating neural stem cell that has been cultured in: a standard cell culture flask such as a T-175 flask; or in a multi-compartment bioreactor for 4 weeks or less.

[0050] A neural stem cell microparticle is a microparticle that is produced by a neural stem cell. Typically, the microparticle is secreted by the neural stem cell. The microparticle of the invention is an exosome. Microparticles from other cells, such as mesenchymal stem cells, are known in the art.

[0051] A "microparticle" is an extracellular vesicle of 30 to 1000 nm diameter that is released from a cell. It is limited by a lipid bilayer that encloses biological molecules. The term "microparticle" is known in the art and encompasses a

number of different species of microparticle, including a membrane particle, membrane vesicle, microvesicle, exosome-like vesicle, exosome, ectosome-like vesicle, ectosome or exovesicle. The invention provides a neural stem cell-derived exosome. The different types of microparticle are distinguished based on diameter, subcellular origin, their density in sucrose, shape, sedimentation rate, lipid composition, protein markers and mode of secretion (i.e. following a signal (inducible) or spontaneously (constitutive)). Four of the common microparticles and their typical distinguishing features are described in Table 1, below.

Table 1: Various Microparticles

Microparticle	Size	Shape	Markers	Lipids	Origin
Microvesicles	100-1000nm	Irregular	Integrins, selectins, CD40 ligand	Phosphatidylserine	Plasma membrane
Exosome-like vesicles	20-50nm	Irregular	TNFR1	No lipid rafts	MVB from other organelles
Exosomes	30-100nm; (<200nm)	Cup shaped	Tetraspanins (e.g. CD63, CD9), Alix , TSG101, ESCRT	Cholesterol, sphingomyelin, ceramide, lipid rafts, phosphatidylserine	Multivesicular endosomes
Membrane particles	50-80nm	Round	CD133, no CD63	Unknown	Plasma membrane

[0052] Microparticles are thought to play a role in intercellular communication by acting as vehicles between a donor and recipient cell through direct and indirect mechanisms. Direct mechanisms include the uptake of the microparticle and its donor cell-derived components (such as proteins, lipids or nucleic acids) by the recipient cell, the components having a biological activity in the recipient cell. Indirect mechanisms include microvesicle-recipient cell surface interaction, and causing modulation of intracellular signalling of the recipient cell. Hence, microparticles may mediate the acquisition of one or more donor cell-derived properties by the recipient cell. It has been observed that, despite the efficacy of stem cell therapies in animal models, the stem cells do not appear to engraft into the host. Accordingly, the mechanism by which stem cell therapies are effective is not clear. Without wishing to be bound by theory, the inventors believe that the microparticles secreted by neural stem cells play a role in the therapeutic utility of these cells and are therefore therapeutically useful themselves.

[0053] The microparticles and stem cells of the disclosure are isolated. The term "isolated" indicates that the microparticle, microparticle population, cell or cell population to which it refers is not within its natural environment. The microparticle, microparticle population, cell or cell population has been substantially separated from surrounding tissue. In some embodiments, the microparticle, microparticle population, cell or cell population is substantially separated from surrounding tissue if the sample contains at least about 75%, in some embodiments at least about 85%, in some embodiments at least about 90%, and in some embodiments at least about 95% microparticles and/or stem cells. In other words, the sample is substantially separated from the surrounding tissue if the sample contains less than about 25%, in some embodiments less than about 15%, and in some embodiments less than about 5% of materials other than the microparticles and/or stem cells. Such percentage values refer to percentage by weight. The term encompasses cells or microparticles which have been removed from the organism from which they originated, and exist in culture. The term also encompasses cells or microparticles which have been removed from the organism from which they originated, and subsequently re-inserted into an organism. The organism which contains the re-inserted cells may be the same organism from which the cells were removed, or it may be a different organism.

[0054] Neural stem cells naturally produce microparticles by a variety of mechanisms, including budding of the plasma membrane (to form membrane vesicles and microvesicles) and as a result of the fusion of intracellular multivesicular bodies (which contain microparticles) with the cell membrane and the release of the microparticles into the extracellular compartment (to secrete exosomes and exosome-like vesicles).

[0055] The neural stem cell that produces the microparticles of the disclosure can be a fetal, an embryonic, or an adult neural stem cell, such as has been described in US5851832, US6777233, US6468794, US5753506 and WO-A-2005121318. The fetal tissue may be human fetal cortex tissue. The cells can be selected as neural stem cells from the differentiation of induced pluripotent stem (iPS) cells, as has been described by Yuan et al. (2011) or a directly induced neural stem cell produced from somatic cells such as fibroblasts (for example by constitutively inducing Sox2, Klf4, and c-Myc while strictly limiting Oct4 activity to the initial phase of reprogramming as recently by Their *et al*, 2012). Human embryonic stem cells may be obtained by methods that preserve the viability of the donor embryo, as is known in the

art (e.g. Chung *et al.* 2008). Such non-destructive methods of obtaining human embryonic stem cell may be used to provide embryonic stem cells from which microparticles of the disclosure can be obtained. Alternatively, microparticles of the disclosure can be obtained from adult stem cells, iPS cells or directly-induced neural stem cells. Accordingly, microparticles of the disclosure can be produced by multiple methods that do not require the destruction of a human embryo or the use of a human embryo as a base material.

[0056] Typically, the neural stem cell population from which the microparticles are produced, is substantially pure. The term "substantially pure" as used herein, refers to a population of stem cells that is at least about 75%, in some embodiments at least about 85%, in some embodiments at least about 90%, and in some embodiments at least about 95% pure, with respect to other cells that make up a total cell population. For example, with respect to neural stem cell populations, this term means that there are at least about 75%, in some embodiments at least about 85%, in some embodiments at least about 90%, and in some embodiments at least about 95% pure, neural stem cells compared to other cells that make up a total cell population. In other words, the term "substantially pure" refers to a population of stem cells of the present disclosure that contain fewer than about 25%, in some embodiments fewer than about 15%, and in some embodiments fewer than about 5%, of lineage committed cells in the original unamplified and isolated population prior to subsequent culturing and amplification.

[0057] A neural stem cell microparticle comprises at least one lipid bilayer which typically encloses a milieu comprising lipids, proteins and nucleic acids. The nucleic acids may be deoxyribonucleic acid (DNA) and/or ribonucleic acid (RNA). RNA may be messenger RNA (mRNA), micro RNA (miRNA) or any miRNA precursors, such as pri-miRNA, pre-miRNA, and/or small nuclear RNA (snRNA).

[0058] A neural stem cell microparticle retains at least one biological function of the stem cell from which it is derived. Biological functions that may be retained include the ability to inhibit cell migration, for example of fibroblasts or fibroblast-like cells, or of a tumour cell such as a glioblastoma cell. In one embodiment, the at least one biological function is that of a neural stem cell that has been cultured in a multi-compartment bioreactor, for at least 10 weeks and optionally no more than 20 weeks. Alternatively the at least one biological function may be that of a neural stem cell-conditioned medium from a neural stem cell that has been cultured in a multi-compartment bioreactor, for at least 10 weeks and optionally no more than 20 weeks. In another embodiment, the at least one biological function is that of a neural stem cell that has been cultured in a T-175 flask under standard conditions.

[0059] Figures 1 and 2 (Example 1) demonstrate that exosomes isolated from the conditioned medium of CTX0E03 cells that have been cultured for 11 weeks have the ability to inhibit fibroblast migration in a transmembrane assay model of cell migration. Accordingly, one biological function that microparticles of the disclosure may retain is the ability to inhibit migration of fibroblast or fibroblast-like cells, for example of normal human dermal fibroblasts (NHDF).

[0060] In contrast, exosomes isolated from the conditioned medium of CTX0E03 cells that have been cultured for 0-6 weeks *promote* cell migration as determined using a scratch/wound closure assay. Examples 1 and 2, Table 2 and Figures 1-3 demonstrate that exosomes isolated from the conditioned medium of CTX0E03 cells that have been cultured for 0-6 weeks retain the ability to close a wound in a "scratch" model of wound healing. The results in Figure 3A show that the migration activity of normal human dermal fibroblasts (NHDF) cultured in CTX0E03 conditioned media is almost the same as the migration activity observed on the addition of purified exosomes.

[0061] The Examples also demonstrate that exosomes isolated from the conditioned medium of CTX0E03 cells that have been cultured for 11 weeks have the ability to promote the destruction of cancer cells by the immune system. Accordingly, one biological function that exosomes disclosed herein may retain is the ability to promote the destruction of cancer cells by the immune system.

[0062] The Examples further demonstrate that exosomes isolated from the conditioned medium of proliferating CTX0E03 cells have the ability to inhibit tumour cell migration. Accordingly, one biological function that exosomes of the invention may retain is the ability to inhibit tumour cell migration of glioblastoma cells.

[0063] Yet further, the Examples demonstrate that exosomes isolated from the conditioned medium of proliferating CTX0E03 cells have the ability to induce differentiation of tumour cells. Accordingly, one biological function that exosomes of the invention may retain is the ability to induce differentiation of glioblastoma cells. Differentiation may readily be determined by known assays, including cell morphology and the presence of stem cell markers (e.g. nestin) and/or differentiated cell markers (e.g. DCX, GFAP). In one embodiment, differentiation is determined by assaying for the stem cell marker nestin (as demonstrated in the Examples). A reduction in nestin expression indicates differentiation of the cell.

Inhibition of cell migration

[0064] Microparticles of the disclosure are able to inhibit cell migration. Typically, the migration of fibroblasts or glioblastoma cells is inhibited. Cell migration assays are known in the art. Two exemplary assays are described below, and are used in the Examples.

[0065] Transmembrane assays (sometimes referred to as "transwell" assays) are known in the art and an exemplary assay is described in Example 1. The assay uses a chamber separated into two compartments by a porous filter

membrane. The cells are seeded on one side of the membrane, while medium containing the purified microparticles is placed on the opposing (lower) side. Example 1 uses fibroblasts, but other cells may be used. After an incubation period (e.g. 6-24 hours), the membrane is fixed and stained to reveal migrated cells (e.g. cell nuclei). The number of cells which have migrated through the pores of the membrane is counted microscopically.

[0066] An alternative transmembrane assay was used in Example 6, as shown in Figure 22. Here, glioblastoma cells are seeded on one side of the porous filter membrane. These cells are either seeded together with 20 μ g/ml microparticles (Figures 22A and C) or have been pre-treated with 10 μ g/ml microparticles for 24 hours (Figure 22B). Medium containing a chemoattractant is placed on the opposing (lower) side. In Example 6, the chemoattractant is FBS. After an incubation period (e.g. 6-24 hours), the membrane is fixed and stained to reveal migrated cells (e.g. cell nuclei). The number of cells which have migrated through the pores of the membrane is again counted microscopically.

[0067] Cell migration is calculated as the number of cells that have migrated through the pores of the membrane in relation to basal conditions (without the microparticles). Inhibition of cell migration in this assay may typically be defined as a decrease in the number of cells that have migrated through the membrane, typically the number of migrated cells is less than 90%, more typically less than 80%, more typically less than 75%, less than 60% or less than 50% of the number of cells that have migrated through the membrane under basal conditions (without the microparticles) after the same incubation period (e.g. 6 hours or 24 hours). As a guideline, inhibition of cell migration is achieved if after 24 hours incubation in the transmembrane assay, the number of human dermal fibroblasts or glioblastoma cells that have migrated through the membrane in the presence of 20 μ g/ml of the microparticles is less than 80% of the number of fibroblasts that migrated under basal conditions (i.e. in the absence of the microparticles).

[0068] In one embodiment, "inhibition of cell migration" is a statistically significant reduction in cell migration of human dermal fibroblasts or glioblastoma cells in a transmembrane assay with a p value of $p < 0.05$, typically $p < 0.001$, in the presence of the microparticles, compared to the migration in the absence of the microparticles. Typically, this is determined after a 24 hour assay incubation period.

[0069] Cell migration may also be determined using an *in vitro* scratch (wound closure) assay, for example the assay of Example 2. Scratch assays were first used as models of wound healing for epithelial or mesenchymal cells. In this assay, cells are seeded into an assay plate and allowed to attach, spread, and form a confluent monolayer. Example 2 uses fibroblasts, but other cells may be used. A pin or needle is used to scratch and remove cells from a discrete area of the confluent monolayer to form a cell-free zone into which cells at the edges of the wound can migrate. Alternatively, a removable insert having a defined shape is placed on contact with the well bottom before the cells are seeded and allowed to form a confluent monolayer excluding the area covered by the insert. The insert is then removed, allowing the cells to migrate onto the newly revealed surface. Using either setup, molecules of interest as potential therapeutics (e.g. the purified microparticles of the disclosure) are added to the well and images of cell movement are captured at regular intervals, for example within a 24-72 hour period, for data analysis.

[0070] Cell migration/wound closure is calculated as the area covered by cells in relation to the initial wound area as determined at 0 hours. Inhibition of cell migration in this assay is typically defined as a decrease in wound closure, typically a wound closure less than 90%, more typically less than 80%, more typically less than 75%, less than 60% or less than 50% of the wound closure observed under basal conditions (without the microparticles) after 24 hours. After 48 hours, the wound closure is typically less than 90% or less than 80% of the wound closure observed using in the absence of the microparticles.

[0071] Inhibition of cell migration may also be defined as delaying a wound closure of 100%, as determined by the scratch assay, by at least 24 hours compared to the wound closure observed under basal conditions. Typically, this delay is achieved by using 2 μ g/ml of the isolated microparticles, as used in Example 2.

[0072] The proteomic analysis in Example 18 indicates that neural stem cell exosomes comprise biological functions associated with the production, packaging, function and degradation of genetic material. Accordingly, in one embodiment, exosomes of the invention retain these functions, typically one or more of RNA polymerase function, RNA degradation function, ribosome function and spliceosome function.

[0073] The microparticle obtained from the neural stem cell has a diameter of 1000nm or less. Typically, the microparticle of the disclosure will have a diameter of 200nm or less, for example 100nm or less. As noted in Table 1 above, microvesicles have a diameter of 100nm to 1000nm. Exosomes are typically defined as having a diameter of 30-100nm, but more recent studies confirm that exosomes can also have a diameter between 100nm and 200nm, (e.g. Katsuda et al, Proteomics 2013 and Katsuda et al, Scientific Reports 2013). Accordingly, exosomes typically have a diameter between 30nm and 150nm. Membrane particles have a diameter of 50nm to 80nm and exosome-like particles have a diameter of 20nm-50nm. The diameter can be determined by any suitable technique, for example electron microscopy or dynamic light scattering. The term microparticle includes, but is not limited to: membrane particle, membrane vesicle, microvesicle, exosome-like vesicle, exosome, ectosome-like vesicle, ectosome or exovesicle.

[0074] Figure 4 panels A-E show the presence in neural stem cells of multivesicular bodies (MVBs) containing exosomes between 30-50nm in diameter, while panel F shows microvesicles >100nm in diameter. Table 21 and Figure 7 (below) show that typical neural stem cell exosomes were measured to have a diameter ranging from approximately

70nm to approximately 150nm, which is consistent with the size of exosomes (from mesenchymal stem cells) described in the art. Accordingly, exosomes of the invention typically have a diameter between 30nm and 200nm, more typically between 50nm and 150nm. As noted above, exosomes are typically positive for the Alix marker (UNIPROT Accession No. Q8WUM4).

[0075] Figure 4F and Table 21 shows the observed size of typical neural stem cell microvesicles, with a mode diameter of approximately 150nm - 200nm, or a median diameter of approximately 180nm - 350nm. Accordingly, microvesicles of the disclosure typically have a diameter between 100 and 1000nm, more typically between 150nm and 350nm.

[0076] Some microparticles of the disclosure express the CD133 surface marker. Other microparticles of the disclosure do not express the CD133 surface marker.

[0077] "Marker" refers to a biological molecule whose presence, concentration, activity, or phosphorylation state may be detected and used to identify the phenotype of a cell.

[0078] Exosomes are endosome-derived lipid microparticles of typically 30-100nm diameter and sometimes between 100nm and 200nm diameter, that are released from the cell by exocytosis. Exosome release occurs constitutively or upon induction, in a regulated and functionally relevant manner. During their biogenesis, exosomes incorporate a wide range of cytosolic proteins (including chaperone proteins, integrins, cytoskeletal proteins and the tetraspanins) and genetic material. Consequently, exosomes are considered to be inter-cellular communication devices for the transfer of proteins, lipids and genetic material between cells, in the parent cell microenvironment and over considerable distance. Although the invention is not bound by this theory, it is possible that the exosomes are responsible for the efficacy of the neural stem cells. Therefore, exosomes from neural stem cells are themselves expected to be therapeutically efficacious.

Microparticles designed to have desired functions

[0079] Microparticles retain at least some of the functions of the stem cells that produce them. Therefore, it is possible to design microparticles by manipulating the stem cell (which can be any stem cell type and is not limited to neural stem cells, although the neural stem cell microparticles of the disclosure are expressly included as an embodiment) to possess one or more desired functions, typically protein or miRNA. The manipulation will typically be genetic engineering, to introduce one or more exogenous coding, non-coding or regulatory nucleic acid sequences into the stem cell. For example, if an exosome containing VEGF and/or bFGF is desired, then the exosome-producing stem cell can be transformed or transfected to express (high levels of) VEGF and/or bFGF, which would then be incorporated into the microparticles produced by that stem cell. Similarly, iPS cells can be used to produce microparticles, and these cells can be designed to produce the proteins and nucleic acids (e.g. miRNA) that are required in the microparticles produced by the iPS cells. The disclosure therefore provides *ad hoc* microparticles, from any stem cell type, that contain a function that is not naturally present in the stem cell from which is produced, i.e. the microparticles (e.g. exosomes) contain one or more exogenous protein or nucleic acid sequences, are not naturally-occurring and are engineered.

[0080] As disclosed herein, isolated or purified microparticles from the conditioned medium of neural stem cells that have been cultured for more than 10 weeks, for example for 11 weeks, and optionally no longer than 20 weeks, are loaded with one or more exogenous nucleic acids, lipids, proteins, drugs or prodrugs which are intended to perform a desired function in a target cell. This does not require manipulation of the stem cell and the exogenous material can optionally be directly added to the microparticles. For example, exogenous nucleic acids can be introduced into the microparticles by electroporation. The microparticles can then be used as vehicles or carriers for the exogenous material. As disclosed herein, microparticles that have been isolated from the cells that produced them may be loaded with exogenous siRNA, typically by electroporation, to produce microparticles that can be deployed to silence one or more pathological genes. In this way, microparticles can be used as vehicles to deliver one or more agents, typically therapeutic or diagnostic agents, to a target cell, for example to enhance or complement their endogenous inhibition of cell migration. An example of this is a neural stem cell exosome comprising exogenous siRNA capable of silencing one or more pathological genes.

Microparticle Marker

[0081] The disclosure provides a population of isolated neural stem cell microparticles, wherein the population essentially comprises only microparticles of the disclosure, i.e. the microparticle population is pure. In many aspects, the microparticle population comprises at least about 80% (in other aspects at least 85%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99%, 99.5%, 99.9% or 100%) of the microparticles of the disclosure.

[0082] The isolated neural stem cell microparticle of the disclosure is characterised in that it has a distinctive expression profile for certain markers and is distinguished from microparticles from other cell types. When a marker is described herein, its presence or absence may be used to distinguish the microparticle. For example, the term "may comprise" or "may express" also discloses the contrary embodiment wherein that marker is not present, e.g. the phrase "the micro-

particle may comprise one or more tetraspanins, typically CD63, CD81, CD9, CD53, CD82 and/or CD37" also describes the contrary embodiment wherein the microparticle may not comprise one or more tetraspanins, typically CD63, CD81, CD9, CD53, CD82 and/or CD37.

[0083] The neural stem cell microparticle of the disclosure is typically considered to carry a marker if at least about 70% of the microparticles of the population, e.g. 70% of the membrane particles, membrane vesicles, microvesicles, exosome-like vesicles, exosomes, ectosome-like vesicles, ectosomes or exovesicles show a detectable level of the marker. In other aspects, at least about 80%, at least about 90% or at least about 95% or at least about 97% or at least about 98% or more of the population show a detectable level of the marker. In certain aspects, at least about 99% or 100% of the population show detectable level of the markers. Quantification of the marker may be detected through the use of a quantitative RT-PCR (qRT-PCR) or through fluorescence activated cell sorting (FACS). It should be appreciated that this list is provided by way of example only, and is not intended to be limiting. Typically, a neural stem cell microparticle of the disclosure is considered to carry a marker if at least about 90% of the microparticles of the population show a detectable level of the marker as detected by FACS.

[0084] The markers described herein are considered to be expressed by a cell of the population of the disclosure, if its expression level, measured by qRT-PCR has a crossing point (Cp) value below or equal to 35 (standard cut off on a qRT-PCR array). The Cp represents the point where the amplification curve crosses the detection threshold, and can also be reported as crossing threshold (ct).

[0085] In one embodiment, the disclosure relates to microparticles produced by a neural stem cell population characterised in that the cells of the population express one or more of the markers Nestin, Sox2, GFAP, β III tubulin, DCX, GALC, TUBB3, GDNF and IDO. In another aspect of the disclosure, the microparticle is an exosome and the population of exosomes expresses one or more of DCX (doublecortin - an early neuronal marker), GFAP (Glial fibrillary acidic protein - an astrocyte marker), GALC, TUBB3, GDNF and IDO.

[0086] The neural stem cell microparticles of the disclosure may express one or more protein markers at a level which is lower or higher than the level of expression of that marker in a mesenchymal stem cell microparticle of the same species. Protein markers that are expressed by the CTX0E03 cell microparticles are identified herein and below. In some embodiments, the microparticles may express a protein marker at a level relative to α tubulin or other such control protein(s). In some embodiments, the microparticles of the disclosure may express that protein at a level of at least ± 1.2 fold change relative to the control protein, typically at least ± 1.5 fold change relative to the control protein, at least ± 2 fold change relative to the control protein or at least ± 3 fold change relative to the control protein. In some embodiments, the microparticles may express a protein marker at a level of between 10^{-2} and 10^{-6} copies per cell relative to α tubulin or other control protein. In some embodiments, the microparticles of the disclosure may express that protein at a level of between 10^{-2} and 10^{-3} copies per cell relative to α tubulin or other control protein.

[0087] The neural stem cell microparticles of the disclosure may express one or more miRNAs (including miRNA precursors) at a level which is lower or higher than the level of expression of that miRNA (including miRNA precursors) in a mesenchymal stem cell microparticle of the same species. miRNA markers that are expressed by the CTX0E03 cell microparticles are identified below. In some embodiments, the microparticles of the disclosure may express the marker miRNA at a level of least ± 1.5 fold change, typically at least ± 2 fold change or at least ± 3 fold change (calculated according to the $\Delta\Delta$ ct method, which is well-known) relative to U6B or 15a, or any other miRNA reference gene, also referred to as an internal control gene.

[0088] The neural stem cell microparticles of the disclosure may express one or more mRNAs at a level which is lower or higher than the level of expression of that mRNA in a mesenchymal stem cell microparticle of the same species. In some embodiments, the microparticles of the disclosure may express the marker mRNA at a level of least ± 1.5 fold change, typically at least ± 2 fold change or at least ± 3 fold change (calculated according to the $\Delta\Delta$ ct method) relative to ATP5B or YWHAZ, or any other reference gene, also referred to as an internal control gene.

[0089] Exosomes of the invention typically express specific integrins, tetraspanins, MHC Class I and/or Class II antigens, CD antigens and cell-adhesion molecules on their surfaces, which may facilitate their uptake by specific cell types. Exosomes contain a variety of cytoskeletal proteins, GTPases, clathrin, chaperones, and metabolic enzymes (but mitochondrial, lysosomal and ER proteins are excluded, so the overall profile does not resemble the cytoplasm). They also contain mRNA splicing and translation factors. Finally, exosomes generally contain several proteins such as HSP70, HSP90, and annexins that are known to play signalling roles yet are not secreted by classical (ER-Golgi) mechanisms.

[0090] The lipid bilayer of an exosome is typically enriched with cholesterol, sphingomyelin and ceramide. Exosomes also express one or more tetraspanin marker proteins. Tetraspanins include CD81, CD63, CD9, CD53, CD82 and CD37. Exosomes can also include growth factors, cytokines and RNA, in particular miRNA. Exosomes typically express one or more of the markers TSG101, Alix, CD109, thy-1 and CD133. Alix (Uniprot accession No. Q8WUM4), TSG101 (Uniprot accession No. Q99816) and the tetraspanin proteins CD81 (Uniprot accession No. P60033) and CD9 (Uniprot accession No. P21926) are characteristic exosome markers.

[0091] Alix is an endosomal pathway marker. Exosomes are endosomal-derived and, accordingly, a microparticle positive for this marker is characterised as an exosome. Exosomes of the invention are typically positive for Alix. Micro-

vesicles of the disclosure are typically negative for Alix.

Microparticle proteome

[0092] Tables 19 and 21 list all proteins detected by mass spectrometry in exosomes and microvesicles, respectively, isolated from CTX0E03 cells cultured for two weeks in an Integra Celine multi-compartment bioreactor. Exosomes of the invention may contain at least a proportion of the proteins identified in Tables 19. Thus, in one embodiment, exosomes of the invention comprise at least 70%, at least 80%, at least 90%, at least 95%, at least 99% or at least 99.5% of the proteins listed in Table 19. Similarly, microvesicles of the disclosure typically comprise at least 70% at least 80%, at least 90%, at least 95%, at least 99% or at least 99.5% of the proteins listed in Table 21. In a further embodiment, the proteome of a microvesicle of the disclosure or exosome of the invention is least 70%, at least 80%, at least 90%, at least 95%, at least 99% or at least 99.5% identical to the proteome provided in Table 19 (exosome) or Table 21 (microvesicle). When determining the protein content of a microparticle or exosome, mass spectrometry is typically used, for example the LC/MS/MS method described in Example 18.

[0093] Tables 20 and 22 show the 100 most abundant proteins detected by mass spectrometry in exosomes and microvesicles, respectively, isolated from CTX0E03 cells cultured for two weeks in an Integra Celine multi-compartment bioreactor. Exosomes of the invention may contain at least a proportion of the proteins identified in Tables 20. Typically, an exosome of the invention comprises the first ten proteins listed in Table 20, more typically the first 20, the first 30, the first 40 or the first 50 proteins listed in Table 20. Similarly, a microparticle of the disclosure typically comprises the first ten proteins listed in Table 22, more typically the first 20, the first 30, the first 40 or the first 50 proteins listed in Table 22. In one embodiment, an exosome of the invention comprises all 100 proteins listed in Table 20. In one embodiment, a microvesicle of the disclosure comprises all 100 proteins listed in Table 22. Typically, the 100 most abundant proteins in an exosome of the invention or microvesicle of the disclosure contain at least 70 of the proteins identified in Table 20 (exosome) or Table 22 (microparticle). More typically, the 100 most abundant proteins in an exosome of the invention or microvesicle of the disclosure contain at least 80, at least 90, at least 95, 96, 97, 98 or 99, or all 100 of the proteins identified in Table 20 (exosome) or Table 22 (microparticle).

Microparticle miRNA content

[0094] Example 17A-C (and the related Figure 13A&B) shows the results of deep sequencing of miRNA present in CTX0E03 cells (standard culture) and in microvesicles and exosomes produced by these cells. This Example shows that, surprisingly, the number of different miRNA species present in the microparticles is greatly reduced compared to the number of different miRNA species present in the cells; the microparticles contain fewer than 120 different miRNAs whereas the cells contain between 450 and 700 miRNA species. The microparticles contain a majority of hsa-miR-1246.

[0095] The data in Example 17 (Tables 5-10) also show that the microparticles are characterised by four main miRNA species, namely hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532. These four miRNAs are the only miRNAs present at a read count of greater than 1000 in the microparticles; these four miRNAs are present in massive excess compared to the other miRNAs in the microparticles. This is in contrast to the profile in the cells, which contain a much greater number of miRNAs present at high (read count greater than 1000) or very high (read count greater than 10,000) levels. Although not bound by theory, the inventors propose that hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532 are selectively trafficked (or otherwise incorporated) into the microparticles and are thought to play a role in the function of the microparticles. A composition may comprise two, three or all four of hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532. This composition is optionally a pharmaceutical composition, comprising a pharmaceutically-acceptable carrier, diluent, vehicle and/or excipient. The pharmaceutical composition is suitable for use in therapy, typically in the same therapies as the microparticles of the disclosure, as noted above.

[0096] Exosomes of the invention may contain at least a proportion of the miRNA species identified in Tables 7-10. In one embodiment, exosomes of the invention comprise at least 70%, at least 80%, at least 90%, at least 95%, at least 99% or at least 99.5% of the miRNAs listed in Tables 8 and 10. Similarly, microvesicles of the disclosure typically comprise at least 70% at least 80%, at least 90%, at least 95%, at least 99% or at least 99.5% of the miRNAs listed in Tables 7 and 9. In a further embodiment, the total miRNA profile of a microvesicle of the disclosure or exosome of the invention is least 70%, at least 80%, at least 90%, at least 95%, at least 99% or at least 99.5% identical to the total miRNA profile provided in Tables 8 and 10 (exosome) or Tables 7 and 9 (microvesicle). When determining the total miRNA profile of a microparticle or exosome, deep sequencing is typically used, for example the method described in Example 17.

[0097] Typically, in one embodiment exosomes, of the invention contain one, two, three or all four of hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532. Each of these miRNA markers is typically present at a read count (optionally determined using the deep sequence technique described in Example 17) of at least 1000 per microparticle. hsa-miR-1246 may optionally have a read count of at least 2000, 5000, 10,000, 20,000, or 25,000 per microparticle.

Hsa-miR-4492 may optionally have a read count of at least 2000, 3000, 4000 or 5000 per microparticle. Hsa-miR-4532 may optionally have a read count of at least 2000 or 3000 per microparticle.

[0098] In one embodiment, each of hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and/or hsa-miR-4532 is present in the microparticle, e.g. exosome, at a higher read count than is present in the cell that produced the microparticle. In particular, miR-1246 typically has a read count in the microparticle at least twice the read count in the cell, more typically at least 4, 5, 6, 7, or 8 times the read count in the cell, and optionally 10, 15 or 20 times the read count in the cell.

[0099] In one embodiment, exosomes of the invention contain hsa-let-7a-5p, hsa-miR-92b-3p, hsa-miR-21-5p, hsa-miR-92a-3p, hsa-miR-10a-5p, hsa-100-5p and/or hsa-99b-5p at a lower read count than is present in the cell that produced the microparticle. Typically, each of these miRNAs has a read count of less than 1000 in the microparticles of the invention, more typically less than 100, for example less than 50. Optionally, exosomes of the invention contain hsa-let-7a-5p at a read count of less than 50 or less than 25.

[0100] In one embodiment, exosomes of the invention contain fewer than 150 types of miRNA (i.e. different miRNA species) when analysed by deep sequencing, typically fewer than 120 types of miRNA.

[0101] In one embodiment, hsa-miR-1246 is the most abundant miRNA in the exosomes of the invention (optionally determined using the deep sequence technique described in Example 17). Typically, at least 40% of the total count of miRNA in exosomes of the invention is hsa-miR-1246. Typically, at least 50% of the total count of miRNA in exosomes of the invention is hsa-miR-1246.

[0102] hsa-miR-4492 is typically the second-most abundant miRNA in the exosomes of the invention. Typically, at least 3% of the total count of miRNA in exosomes of the invention is hsa-miR-4492. More typically, at least 4% of the total count of miRNA in exosomes of the invention is hsa-miR-4492.

[0103] Typically, at least 2% of the total count of miRNA in exosomes of the invention is hsa-miR-4532.

[0104] Typically, at least 1% of the total count of miRNA in exosomes of the invention is hsa-miR-4488.

[0105] In one embodiment exosomes of the invention contain one or both of hsa-miR-4508, hsa-miR-4516 at a level at least 0.1% of the total miRNA content of the particle.

[0106] One or more of hsa-miR-3676-5p, hsa-miR-4485, hsa-miR-4497, hsa-miR-21-5p, hsa-miR-3195, hsa-miR-3648, hsa-miR-663b, hsa-miR-3656, hsa-miR-3687, hsa-miR-4466, hsa-miR-4792, hsa-miR-99b-5p and hsa-miR-1973 may be present in the exosomes of the invention.

[0107] Typically, each of hsa-let-7a-5p and hsa-100-5p is present at less than 1%, more typically less than 0.1% or less than 0.05% of the total miRNA count in exosomes of the invention.

[0108] In a typical exosome of the invention, at least 50% of the total count of miRNA is hsa-miR-1246, and less than 0.1% of the total miRNA count is hsa-let-7a-5p.

[0109] In one embodiment, at least 90% of the total count of miRNA in exosomes of the invention comprises hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532. Typically, at least 95% or 96% of the total count of miRNA in exosomes of the invention comprises hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532. Less than 10% of the total miRNA content of these exosomes is an miRNA that is not hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532.

[0110] Combinations of the miRNA embodiments discussed above are provided. For example, an exosome of the invention typically contains each of hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532 at a read count of at least 1000 and contains each of hsa-let-7a-5p, hsa-miR-92b-3p, hsa-miR-21-5p, hsa-miR-92a-3p, hsa-miR-10a-5p, hsa-100-5p and hsa-99b-5p at a read count of less than 100. Typically, at least 90% or at least 95% of the total miRNA in these exosomes is hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532.

[0111] An exosome of the invention typically has hsa-miR-1246 as the most abundant miRNA and hsa-miR-4492 is the second-most abundant miRNA. In this embodiment, at least 40% of the total count of miRNA in exosomes of the invention is hsa-miR-1246 and at least 3% of the total count of miRNA in the microparticle is hsa-miR-4492. At least 2% of the total count of miRNA in these exosomes is hsa-miR-4532 and at least 1% of the total count of miRNA in these exosomes is hsa-miR-4488. Each of hsa-let-7a-5p and hsa-100-5p is present at less than 0.1% of the total miRNA count in these exosomes.

[0112] Plotting the deep sequencing results in the exosomes and microvesicles as relative fold change compared to the cells confirms that hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532 are significantly upregulated in the exosomes and microvesicles compared to the cells.

[0113] This comparison also shows that miRNA hsa-miR-3195 is the miRNA that is most upregulated, in both exosomes and microvesicles. Although the absolute reads of hsa-miR-3195 are in the range of -40 for exosomes and microvesicles, there is no hsa-miR-3195 detected in the cells. Accordingly, hsa-miR-3195 is uniquely found in the exosomes and microvesicles and, optionally, in one embodiment, an exosome of the invention comprises hsa-miR-3195.

[0114] In one embodiment, microparticles of the disclosure comprise one or more of the following miRNA precursors:

AC079949.1 (SEQ ID NO:738)

GGCCGCGCCCCGTTTCCCAGGACAAAGGGCACTCCGCACCGGACCCTGGTCCCAGCG;

AP000318.1 (SEQ ID NO: 739)
CCCACTCCCTGGCGCCGCTTGTGGAGGGCCCAAGTCTTCTGATTGAGGCCCAACCCGTGGAAG;

5

10

and

15

20

25

30

35

40

45

50

EP 3 033 418 B1

(continued)

5

10

15

20

25

30

35

40

45

50

55

hsa-miR-101
hsa-miR-103a
hsa-miR-106b
hsa-miR-10a
hsa-miR-10b
hsa-miR-124
hsa-miR-125a-5p
hsa-miR-125b
hsa-miR-126
hsa-miR-127-5p
hsa-miR-128
hsa-miR-129-5p
hsa-miR-130a
hsa-miR-132
hsa-miR-134
hsa-miR-137
hsa-miR-141
hsa-miR-146b-5p
hsa-miR-150
hsa-miR-155
hsa-miR-15a
hsa-miR-15b
hsa-miR-16
hsa-miR-17
hsa-miR-181a
hsa-miR-182
hsa-miR-183
hsa-miR-185
hsa-miR-18a
hsa-miR-18b
hsa-miR-192
hsa-miR-194
hsa-miR-195
hsa-miR-196a
hsa-miR-205
hsa-miR-20a
hsa-miR-20b
hsa-miR-21
hsa-miR-210

(continued)

5

10

15

20

25

30

hsa-miR-214
hsa-miR-218
hsa-miR-219-5p
hsa-miR-22
hsa-miR-222
hsa-miR-23b
hsa-miR-24
hsa-miR-26a
hsa-miR-301a
hsa-miR-302a
hsa-miR-302c
hsa-miR-33a
hsa-miR-345
hsa-miR-375
hsa-miR-378
hsa-miR-424
hsa-miR-7
hsa-miR-9
hsa-miR-92a
hsa-miR-93
hsa-miR-96
hsa-miR-99a

35

[0118] In one embodiment, the CTX0E03 microparticles contain 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 15, 20, 25, 30 or more of the following miRNAs (which are selected from the list above):

40

45

50

55

hsa-let-7g
hsa-miR-101
hsa-miR-10a
hsa-miR-10b
hsa-miR-126
hsa-miR-128
hsa-miR-129-5p
hsa-miR-130a
hsa-miR-134
hsa-miR-137
hsa-miR-155
hsa-miR-15a
hsa-miR-15b
hsa-miR-16

(continued)

5

10

15

20

25

30

35

hsa-miR-17
hsa-miR-182
hsa-miR-183
hsa-miR-185
hsa-miR-18b
hsa-miR-192
hsa-miR-194
hsa-miR-195
hsa-miR-20a
hsa-miR-20b
hsa-miR-210
hsa-miR-218
hsa-miR-301a
hsa-miR-302a
hsa-miR-302c
hsa-miR-345
hsa-miR-375
hsa-miR-378
hsa-miR-7
hsa-miR-9
hsa-miR-93
hsa-miR-96
hsa-miR-99a

miRNAs present in exosomes from cells cultured in a bioreactor for longer periods

[0119] Examples 17D and 17E (in particular Figures 13D to 13H and Tables E2 to E4) demonstrate that hsa-miR-1246, hsa-miR-4492, hsa-miR-4532, and hsa-miR-4488 are still present in exosomes isolated from CTX0E03 cells that have been cultured in a bioreactor for six weeks. hsa-miR-4492, hsa-miR-4532, and hsa-miR-4488 are shown to be almost absent in exosomes isolated from CTX0E03 cells that have been cultured in a bioreactor for eleven weeks.

[0120] Exosomes of the invention may contain at least a proportion of the miRNA species identified in Table E3, or at least a proportion of the miRNA species identified in Table E4.

[0121] Hsa-miR-181a-5p, hsa-miR-1246, hsa-miR-127-3p, hsa-miR-21-5p, and hsa-miR-100-5p are shown to be the top 5 miRNAs present in the EXO 6W sample. Accordingly, in one embodiment, exosomes of the invention comprise 1, 2, 3, 4 or 5 of these miRNAs. In another embodiment, exosomes of the invention comprise at least 50%, at least 60%, at least 70%, at least 80%, at least 90%, or all of the miRNAs listed in Table E3. When determining the total miRNA profile of a microparticle or exosome, deep sequencing is typically used, for example the method described in Example 17.

[0122] Hsa-miR-181a-5p, hsa-let-7a-5p, hsa-let-7f-5p, hsa-miR-92b-3p, and hsa-miR-9-5p are shown to be the top 5 miRNAs present in EXO 11W samples. Accordingly, in one embodiment, exosomes of the invention comprise 1, 2, 3, 4 or 5 of these miRNAs. In another embodiment, exosomes of the invention comprise at least 50%, at least 60%, at least 70%, at least 80%, at least 90%, or all of the miRNAs listed in Table E4. When determining the total miRNA profile of a microparticle or exosome, deep sequencing is typically used, for example the method described in Example 17.

[0123] Hsa-miR-486-5p is observed to be shuttled into all three of the samples of exosomes obtained from CTX0E03 cells that have been cultured in a bioreactor for six weeks. Accordingly, in one embodiment, exosomes of the invention comprise hsa-miR-486-5p.

Individual miRNAs are able to reduce cell proliferation and have therapeutic utility

[0124] The data in the Example 20 and Figure 23 show that each of the four main miRNA species identified in neural stem cell microparticles, namely hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532, significantly reduced cell proliferation in glioma proliferation assays. In addition to these data supporting the therapeutic efficacy of the microparticles that contain these miRNAs, these data also show that each of these individual miRNAs is therapeutically useful on *its own*. In one embodiment, the individual miRNA is useful in the treatment of cancer (optionally glioblastoma), as described below.

Proteins detected by a dot-blot

[0125] Example 13 shows proteins present in microparticles produced by the CTX0E03 cells, as detected by a dot-blot. Microparticles of the invention may typically contain 1, 2, 3, 4, 5, 6, 7, 8, 9, 10 or all of the following proteins:

EDA-A2
Galectin-3
IGFBP-2
IGFBP-rp1/IGFBP-7
IL-1a
LECT2
MCP-1
SPARC
TIMP-1
Thrombospondin-1
VEGF

[0126] Galectin-3 and Thrombospondin-1 are also identified as present in exosomes and microvesicles in Example 18. TIMP-1 is identified in Example 18 as being present in exosomes. Microparticles of the invention may contain one or more of Galectin-3, Thrombospondin and TIMP-1.

[0127] Example 13 also shows that the microparticles produced by the CTX0E03 cells may also express 1, 2, 3, 4 or 5 of the following proteins:

EGF-R/ErbB1
MDC
Endostatin
Follistatin
Csk

[0128] EGF-R and Csk are also identified as present in exosomes and microvesicles in Example 18.

Neural Stem cells in multi-compartment bioreactor culture

[0129] As shown in Example 15 and Figure 11 below, after multi-compartment bioreactor culture for three weeks, neural stem cells express a number of markers at significantly higher levels than neural stem cells cultured according to standard procedure in a standard single-compartment T175 flask. Neural stem cells cultured for even longer periods, e.g. at least 10 weeks, may also express a number of these markers at significantly higher levels than neural stem cells cultured according to standard procedure in a standard single-compartment T175 flask or neural stem cells cultured in a multi-compartment bioreactor culture for three weeks. In one embodiment, microparticles of the disclosure are isolated

from NSCs that have been cultured, typically in a multi-compartment bioreactor, for at least 10 weeks, typically at least 11 weeks, at least 12 weeks, at least 13 weeks, at least 14 weeks or at least 15 weeks. Optionally, the NSCs have been cultured for no more than 20 weeks, e.g. between 10 and 20 weeks, between 11 and 20 weeks, between 12 and 20 weeks, between 13 and 20 weeks, between 14 and 20 weeks or between 15 and 20 weeks.

[0130] CTX0E03 neural stem cells cultured for three weeks in a multi-compartment bioreactor express DCX, GALC, GFAP, TUBB3, GDNF and IDO at a higher level than neural stem cells cultured in a standard single-compartment T175 cell culture. Neural stem cells cultured for even longer periods, e.g. at least 10 weeks, may also express a number of these markers at significantly higher levels than neural stem cells cultured according to standard procedure in a standard single-compartment T175 flask or, optionally, than neural stem cells cultured in a multi-compartment bioreactor culture for three weeks. Accordingly neural stem cells that produce microparticles of the disclosure may express one or more of DCX, GALC, GFAP, TUBB3, GDNF and IDO. Cells cultured in a two-compartment bioreactor typically show increased expression of one or more of DCX, GALC, GFAP, TUBB3, GDNF and IDO compared to the stem cells cultured under standard conditions for three weeks. The expression level of these markers in the multi-compartment bioreactor-cultured cells is typically significantly higher than in the cells cultured in a standard single-compartment T175 culture flask. Typically, a stem cell cultured in a multi-compartment bioreactor, that produces microparticles of the disclosure, expresses one or more of DCX1, GALC, GFAP, TUBB3, GDNF or IDO at a level least 2 fold higher than in CTX0E03 cells cultured in a T-175 flask according to standard culture procedure.

[0131] In one embodiment, microparticles, typically exosomes, are obtained from neural stem cells that show increased expression of one or more of DCX, GALC, GFAP, TUBB3, GDNF and IDO compared to the stem cells cultured under standard conditions or, optionally than in a multi-compartment bioreactor culture for three weeks. For example, microparticles can be obtained from freshly filtered conditioned medium collected from Integra CeLLine bioreactor cultured neural stem cells.

[0132] The upregulated markers include DCX (doublecortin - an early neuronal marker), GFAP (Glial fibrillary acidic protein - an astrocyte marker), GALC, TUBB3, GDNF and IDO. CTX0E03 cells are able to differentiate into 3 different cell types: neurons, astrocytes and oligodendrocytes. The high levels of DCX and GFAP after only three weeks in a multi-compartment bioreactor indicates that the cultured stem cells have partially differentiated and have entered the neuronal (DCX+ cells) and/or astrocytic (GFAP+ cells) lineage. Accordingly, in one embodiment the disclosure provides a microparticle that inhibits cell migration, produced by a neural stem cell population that expresses (i) one or more markers associated with a neuronal lineage, typically DCX and/or (ii) one or more markers associated with an astrocytic lineage, typically GFAP. These cells may optionally have been cultured for at least 10 weeks in a multi-compartment bioreactor. In another embodiment, the disclosure provides neural stem cell microparticles, typically exosomes, that express (i) one or more markers associated with a neuronal lineage, typically DCX and/or (ii) one or more markers associated with an astrocytic lineage, typically GFAP. These cells, or the microparticles (typically exosomes) derived from these cells, express DCX and/or GFAP at a higher level than the corresponding stem cells in standard (T-175) culture or, optionally, than the cells cultured in a multi-compartment bioreactor for three weeks. Typically, these cells or microparticles express DCX and/or GFAP at a level at least 2 fold more than the stem cells in standard culture, more typically at least 2.5 fold more than the corresponding stem cells in standard culture (or cultured in a multi-compartment bioreactor culture for three weeks), at least 5 fold more than the corresponding stem cells in standard culture (or cultured in a multi-compartment bioreactor culture for three weeks), at least 7.5 fold more than the corresponding stem cells in standard culture (or cultured in a multi-compartment bioreactor culture for three weeks) or at least 10 fold more than the corresponding stem cells in standard culture (or cultured in a multi-compartment bioreactor culture for three weeks). For expression of DCX, the fold change in the cells or microparticles compared to the corresponding stem cells in standard (T-175) culture (or cultured in a multi-compartment bioreactor culture for three weeks) can optionally be at least 20 fold, at least 50 fold, at least 100 fold, at least 500 fold or at least 1000 fold more than the standard stem cells (or cells cultured in a multi-compartment bioreactor culture for three weeks).

[0133] The term "bioreactor" is to be given its usual meaning in the art, i.e. an apparatus used to carry out a bioprocess. The bioreactors described herein are suitable for use in stem cell culture. Simple bioreactors for cell culture are single compartment flasks, such as the commonly-used T-175 flask (e.g. the BD Falcon™ 175 cm² Cell Culture Flask, 750 ml, tissue-culture treated polystyrene, straight neck, blue plug-seal screw cap, BD product code 353028).

[0134] Bioreactors can have multiple compartments, as is known in the art. These multi-compartment bioreactors typically contain at least two compartments separated by one or more membranes or barriers that separate the compartment containing the cells from one or more compartments containing gas and/or culture medium. Multi-compartment bioreactors are well-known in the art. An example of a multi-compartment bioreactor is the Integra CeLLine bioreactor, which contains a medium compartment and a cell compartment separated by means of a 10 kDa semipermeable membrane; this membrane allows a continuous diffusion of nutrients into the cell compartment with a concurrent removal of any inhibitory waste product. The individual accessibility of the compartments allows to supply cells with fresh medium without mechanically interfering with the culture. A silicone membrane forms the cell compartment base and provides an optimal oxygen supply and control of carbon dioxide levels by providing a short diffusion pathway to the cell com-

partment. Any multi-compartment bioreactor may be used according to the disclosure. As shown in the Examples below, CTX0E03 cells that have been cultured in the Integra CeLLine AD1000 bioreactor for 11 weeks produce microparticles that are able to inhibit cell migration.

[0135] Example 16, Table 4 and Figure 12 show that the miRNA content of exosomes produced by neural stem cells that have been cultured in a multi-compartment bioreactor, for three weeks, is different from the miRNA content of stem cells cultured in standard T-175 flasks and from microparticles produced by the neural stem cells cultured in a single-compartment T175 culture flask for three weeks. The miRNA content of exosomes of the invention may also differ from the miRNA content of stem cells cultured in standard T-175 or microparticles derived therefrom. In one embodiment, the invention provides a microparticle, typically an exosome, wherein at least two, three, four, five, six or seven miRNAs are up or down regulated compared to in the corresponding stem cells cultured in standard T-175 flasks, as calculated by Fold Regulation (see Example 16), and wherein the microparticle inhibits cell migration. The Fold Regulation of each miRNA is optionally at least two-fold up or down.

[0136] In one embodiment, neural stem cell exosomes of the invention express one, two, three, four, five, six or seven of the following miRNAs at a higher level than is expressed in the corresponding stem cells cultured in standard T-175 flasks, as calculated by Fold Regulation (where an asterisk indicates an miRNA where at least a two-fold regulation increase is preferred):

hsa-miR-146b-5p*
hsa-let-7c*
hsa-miR-99a*
hsa-miR-132*
hsa-miR-378*
hsa-miR-181a*
hsa-let-7b*

[0137] In one embodiment, neural stem cell exosomes of the invention express one, two, three, four, five, six, seven, eight, nine, ten or more of the following miRNAs at a lower level than is expressed in the corresponding stem cells cultured in standard T-175 flasks, as calculated by Fold Regulation (where an asterisk indicates an miRNA where at least a two-fold regulation decrease is preferred):

hsa-miR-7*
hsa-miR-106b*
hsa-miR-101*
hsa-miR-302a*
hsa-miR-301a*
hsa-miR-183*
hsa-miR-219-5p*
hsa-miR-18a*
hsa-miR-15a*
hsa-miR-182*
hsa-miR-33a*
hsa-miR-96*
hsa-miR-18b*

[0138] In a further embodiment, NSC exosomes of the invention comprise (i) an increased level of at least one, two, three, four, five, six or seven of the miRNAs indicated above as being increased in exosomes compared to the corresponding cells in standard culture and (ii) a decreased level of at least one, two, three, four, five, six, seven, eight, nine, ten or more or more of the miRNAs indicated above as being decreased in exosomes compared to the corresponding

cells in standard culture. For example, a neural stem cell exosome may contain a fold-regulation increase in three or more or more of the miRNAs indicated above as being increased in exosomes compared to the corresponding cells in standard culture and a fold-regulation decrease in three or more of the miRNAs indicated above as being decreased in exosomes compared to the corresponding cells in standard culture. In another exemplary embodiment, a neural stem cell exosome may contain a fold-regulation increase in five or more of the miRNAs indicated above as being increased in exosomes compared to the corresponding cells in standard culture and a fold-regulation decrease in five or more of the miRNAs indicated above as being decreased in exosomes compared to the corresponding cells in standard culture.

[0139] The term "expressed" is used to describe the presence of a marker within a cell or microparticle. In order to be considered as being expressed, a marker must be present at a detectable level. By "detectable level" is meant that the marker can be detected using one of the standard laboratory methodologies such as qRT-PCR, or qPCR, blotting, Mass Spectrometry or FACS analysis. A gene is considered to be expressed by a cell or microparticle of the population of the invention if expression can be reasonably detected at a crossing point (cp) values below or equal 35. The terms "express" and "expression" have corresponding meanings. At an expression level below this cp value, a marker is considered not to be expressed. The comparison between the expression level of a marker in a stem cell or microparticle of the disclosure, and the expression level of the same marker in another cell or microparticle, such as for example an mesenchymal stem cell, may preferably be conducted by comparing the two cell/microparticle types that have been isolated from the same species. Preferably this species is a mammal, and more preferably this species is human. Such comparison may conveniently be conducted using a reverse transcriptase polymerase chain reaction (RT-PCR) experiment.

[0140] As used herein, the term "significant expression" or its equivalent terms "positive" and "+" when used in regard to a marker shall be taken to mean that, in a cell or microparticle population, more than 20%, preferably more than, 30%, 40%, 50%, 60%, 70%, 80%, 90%, 95%, 98%, 99% or even all of the cells of the cells/microparticles express said marker.

[0141] As used herein, "negative" or "-" as used with respect to markers shall be taken to mean that, in a cell or microparticle population, less than 20%, 10%, preferably less than 9%, 8%, 7%, 6%, 5%, 4%, 3%, 2%, 1 % or none of the cells/microparticles express said marker.

[0142] Expression of microparticle surface markers may be determined, for example, by means of flow cytometry and/or FACS for a specific cell surface marker using conventional methods and apparatus (for example a Beckman Coulter Epics XL FACS system used with commercially available antibodies and standard protocols known in the art) to determine whether the signal for a specific microparticle surface marker is greater than a background signal. The background signal is defined as the signal intensity generated by a non-specific antibody of the same isotype as the specific antibody used to detect each surface marker. For a marker to be considered positive the specific signal observed is typically more than 20%, preferably stronger than 30%, 40%, 50%, 60%, 70%, 80%, 90%, 100%, 500%, 1000%, 5000%, 10000% or above, greater relative to the background signal intensity. Alternative methods for analysing expression of microparticle surface markers of interest include visual analysis by electron microscopy using antibodies against cell-surface markers of interest.

[0143] "Fluorescence activated cell sorting (FACS)" is a method of cell purification based on the use of fluorescent labelled antibodies. The antibodies are directed to a marker on the cell surface, and therefore bind to the cells of interest. The cells are then separated based upon the fluorescent emission peak of the cells.

[0144] Microparticle markers (including surface and intracellular proteins) can also be analysed by various methods known to one skilled in the art to assay protein expression, including but not limited to gel electrophoresis followed by western blotting with suitable antibodies, immunoprecipitation followed by electrophoretic analysis, and/or electron microscopy as described above, with microparticle permeabilisation for intraparticle markers. For example, expression of one or more tetraspanins may be assayed using one or more of the above methods or any other method known to one skilled in the art. RNA levels may also be analysed to assess marker expression, for example qRT-PCR.

Microparticle Function

[0145] As noted above, a neural stem cell microparticle typically retains at least one biological function of the stem cell from which it is derived. Biological functions that may be retained include the ability to: inhibit cell migration, for example of fibroblast or fibroblast-like cells, or of tumour cells such as glioblastoma cells; inhibit wound healing, for example in a scratch assay; or treat a disease or condition that involves or is characterised by undesirable or excessive cell migration, such as cancer, fibrosis, atherosclerosis or rheumatoid arthritis.

[0146] In one embodiment, the at least one biological activity is that of a neural stem cell that has been cultured, typically in a multi-compartment bioreactor, for at least 10 weeks and optionally no more than 20 weeks. Alternatively the at least one biological activity may be that of a neural stem cell-conditioned medium from a neural stem cell that has been cultured, typically in a multi-compartment bioreactor, for at least 10 weeks and optionally no more than 20 weeks. Figures 1 and 2 (Example 1) demonstrate that exosomes isolated from the conditioned medium of CTX0E03 cells that have been cultured in a CeLLine bioreactor for 11 weeks have the ability to inhibit fibroblast migration in a transmembrane assay model of cell migration. Accordingly, one biological function that microparticles of the disclosure may retain is the

ability to inhibit migration of fibroblast or fibroblast-like cells, for example of normal human dermal fibroblasts (NHDF).

[0147] Example 2, Table 2 and Figure 3 demonstrate that CTX0E03 stem cell exosomes, obtained from cells cultured for 2 weeks and 6 weeks, retain the ability to close a wound in a "scratch" model of wound healing. The results show that the migration activity of normal human dermal fibroblasts (NHDF) cultured in CTX0E03 conditioned media is almost the same as the migration activity observed on the addition of purified exosomes. In contrast, microparticles of the disclosure are able to inhibit cell migration. Accordingly, one biological function that microparticles of the disclosure may retain is the ability to inhibit migration activity of normal human dermal fibroblasts (NHDF). NHDF migration assays are known in the art. Stimulation of NHDF migration may be determined using an *in vitro* scratch (wound closure) assay, for example the assay of Example 2. Wound closure is calculated as the area covered by NHDF cells in relation to the initial wound area as determined at 0 hours. Inhibition of NHDF migration in this assay is typically defined as a decrease in wound closure, as defined above.

[0148] CTX0E03 cells are known to inhibit T cell activation in a PBMC assay and, in one embodiment, the microparticles of the disclosure retain this ability to inhibit T cell activation in a PBMC assay. PBMC assays are well-known to the skilled person and kits for performing the assay are commercially available.

[0149] The proteomic analysis in Example 18 indicates that neural stem cell exosomes comprise biological functions associated with the production, packaging, function and degradation of genetic material. Accordingly, in one embodiment, exosomes of the invention retain these functions, typically one or more of RNA polymerase function, RNA degradation function, ribosome function and spliceosome function.

Immunogenicity

[0150] The (allogeneic) neural stem cell microparticles of the disclosure typically either do not trigger an immune response *in vitro* or *in vivo* or trigger an immune response which is substantially weaker than that which would be expected to be triggered upon injection of an allogeneic stem cell population into a patient. In certain aspects of the disclosure, the neural stem cell microparticles are considered not to trigger an immune response if at least about 70% of the microparticles do not trigger an immune response. In some embodiments, at least about 80%, at least about 90% or at least about 95%, 99% or more of the microparticles do not trigger an immune response. Preferably the microparticles of the disclosure do not trigger an antibody mediated immune response or do not trigger a humoral immune response. More preferably the microparticles of the disclosure do not trigger either an antibody mediated response or a humoral immune response *in vitro*. More preferably still, the microparticles of the disclosure do not trigger a mixed lymphocyte immune response. It will be understood by one skilled in the art that the ability of the cells of the disclosure to trigger an immune response can be tested in a variety of ways.

[0151] CTX0E03 cells transplanted in a rodent model of limb ischemia have been previously demonstrated a faster and transient up-regulation of host genes involved in angiogenesis, such as CCL11, CCL2, CXCL1, CXCL5, IGF1, IL1 β , IL6, HGF, HIF1 α , bFGF, VEGFA, and VEGFC, compared to vehicle treated controls. hNSC treatment transiently elevates host innate immune and angiogenic responses and accelerates tissue regeneration.

[0152] The CTX0E03 cell line has been previously demonstrated, using a human PBMC assay, not to be immunogenic. Accordingly, microparticles produced by CTX0E03 cells are also expected to be non-immunogenic. The lack of immunogenicity allows the microparticles to avoid clearance by the host/patient immune system and thereby exert their therapeutic effect without a deleterious immune and inflammatory response.

Neural Stem Cells

[0153] The neural stem cell that produces the microparticle may be a stem cell line, i.e. a culture of stably dividing stem cells. A stem cell line can be grown in large quantities using a single, defined source. Immortalisation may arise from a spontaneous event or may be achieved by introducing exogenous genetic information into the stem cell which encodes immortalisation factors, resulting in unlimited cell growth of the stem cell under suitable culture conditions. Such exogenous genetic factors may include the gene "myc", which encodes the transcription factor Myc. The exogenous genetic information may be introduced into the stem cell through a variety of suitable means, such as transfection or transduction. For transduction, a genetically engineered viral vehicle may be used, such as one derived from retroviruses, for example lentivirus.

[0154] Additional advantages can be gained by using a conditionally immortalised stem cell line, in which the expression of the immortalisation factor can be regulated without adversely affecting the production of therapeutically effective microparticles. This may be achieved by introducing an immortalisation factor which is inactive unless the cell is supplied with an activating agent. Such an immortalisation factor may be a gene such as c-mycER. The c-MycER gene product is a fusion protein comprising a c-Myc variant fused to the ligand-binding domain of a mutant estrogen receptor. C-MycER only drives cell proliferation in the presence of the synthetic steroid 4-hydroxytamoxifen (4-OHT) (Littlewood et al.1995). This approach allows for controlled expansion of neural stem cells *in vitro*, while avoiding undesired *in vivo*

effects on host cell proliferation (e.g. tumour formation) due to the presence of c-Myc or the gene encoding it in micro-particles derived from the neural stem cell line. A suitable c-mycER conditionally immortalized neural stem cell is described in United States Patent 7416888. The use of a conditionally immortalised neural stem cell line therefore provides an improvement over existing stem cell microparticle isolation and production.

[0155] Preferred conditionally-immortalised cell lines include the CTX0E03, STR0C05 and HPC0A07 neural stem cell lines, which have been deposited at the European Collection of Animal Cultures (ECACC), Vaccine Research and Production laboratories, Public Health Laboratory Services, Porton Down, Salisbury, Wiltshire, SP4 0JG, with Accession No. 04091601 (CTX0E03); Accession No.04110301 (STR0C05); and Accession No.04092302 (HPC0A07).

[0156] The derivation and provenance of these cells is described in EP1645626 B1. The advantages of these cells are retained by microparticles produced by these cells.

[0157] The cells of the CTX0E03 cell line may be cultured in the following culture conditions:

- Human Serum Albumin 0.03%
- Transferrin, Human 5µg/ml
- Putrescine Dihydrochloride 16.2 µg/ml
- Insulin Human recombinant 5 µg/ml
- Progesterone 60 ng/ml
- L-Glutamine 2 mM
- Sodium Selenite (selenium) 40 ng/ml

[0158] Plus basic Fibroblast Growth Factor (10 ng/ml), epidermal growth factor (20 ng/ml) and 4-hydroxytamoxifen 100nM for cell expansion. The cells can be differentiated by removal of the 4-hydroxytamoxifen. Typically, the cells can either be cultured at 5% CO₂/37°C or under hypoxic conditions of 5%, 4%, 3%, 2% or 1% O₂. These cell lines do not require serum to be cultured successfully. Serum is required for the successful culture of many cell lines, but contains many contaminants including its own exosomes. A further advantage of the CTX0E03, STR0C05 or HPC0A07 neural stem cell lines, or any other cell line that does not require serum, is that the contamination by serum is avoided.

[0159] The cells of the CTX0E03 cell line (and microparticles derived from these cells) are multipotent cells originally derived from 12 week human fetal cortex. The isolation, manufacture and protocols for the CTX0E03 cell line is described in detail by Sinden, et al. (U.S. Pat. 7,416,888 and EP1645626 B1). The CTX0E03 cells are not "embryonic stem cells", i.e. they are not pluripotent cells derived from the inner cell mass of a blastocyst; isolation of the original cells did not result in the destruction of an embryo. In growth medium CTX0E03 cells are nestin-positive with a low percentage of GFAP positive cells (i.e. the population is negative for GFAP).

[0160] CTX0E03 is a clonal cell line that contains a single copy of the c-mycER transgene that was delivered by retroviral infection and is conditionally regulated by 4-OHT (4-hydroxytamoxifen). The C-mycER transgene expresses a fusion protein that stimulates cell proliferation in the presence of 4-OHT and therefore allows controlled expansion when cultured in the presence of 4-OHT. This cell line is clonal, expands rapidly in culture (doubling time 50-60 hours) and has a normal human karyotype (46 XY). It is genetically stable and can be grown in large numbers. The cells are safe and non-tumorigenic. In the absence of growth factors and 4-OHT, the cells undergo growth arrest and differentiate into neurons and astrocytes. Once implanted into an ischemia-damaged brain, these cells migrate only to areas of tissue damage.

[0161] The development of the CTX0E03 cell line has allowed the scale-up of a consistent product for clinical use. Production of cells from banked materials allows for the generation of cells in quantities for commercial application (Hodges et al, 2007).

[0162] Pollock *et al* 2006 describes that transplantation of CTX0E03 in a rat model of stroke (MCAo) caused statistically significant improvements in both sensorimotor function and gross motor asymmetry at 6-12 weeks post-grafting. These data indicate that CTX0E03 has the appropriate biological and manufacturing characteristics necessary for development as a therapeutic cell line.

[0163] Stevanato et al 2009 confirms that CTX0E03 cells downregulated c-mycERTAM transgene expression both *in vitro* following EGF, bFGF and 4-OHT withdrawal and *in vivo* following implantation in MCAo rat brain. The silencing of the c-mycERTAM transgene in vivo provides an additional safety feature of CTX0E03 cells for potential clinical application.

[0164] Smith *et al* 2012 describe preclinical efficacy testing of CTX0E03 in a rat model of stroke (transient middle cerebral artery occlusion). The results indicate that CTX0E03 implants robustly recover behavioural dysfunction over a 3 month time frame and that this effect is specific to their site of implantation. Lesion topology is potentially an important factor in the recovery, with a stroke confined to the striatum showing a better outcome compared to a larger area of damage.

[0165] Neural retinal stem cell lines (for example as described in US 7514259) may also be used according to the invention.

[0166] The term "culture medium" or "medium" is recognized in the art, and refers generally to any substance or

preparation used for the cultivation of living cells. The term "medium", as used in reference to a cell culture, includes the components of the environment surrounding the cells. Media may be solid, liquid, gaseous or a mixture of phases and materials. Media include liquid growth media as well as liquid media that do not sustain cell growth. Media also include gelatinous media such as agar, agarose, gelatin and collagen matrices. Exemplary gaseous media include the gaseous phase to which cells growing on a petri dish or other solid or semisolid support are exposed. The term "medium" also refers to material that is intended for use in a cell culture, even if it has not yet been contacted with cells. In other words, a nutrient rich liquid prepared for culture is a medium. Similarly, a powder mixture that when mixed with water or other liquid becomes suitable for cell culture may be termed a "powdered medium". "Defined medium" refers to media that are made of chemically defined (usually purified) components. "Defined media" do not contain poorly characterized biological extracts such as yeast extract and beef broth. "Rich medium" includes media that are designed to support growth of most or all viable forms of a particular species. Rich media often include complex biological extracts. A "medium suitable for growth of a high density culture" is any medium that allows a cell culture to reach an OD600 of 3 or greater when other conditions (such as temperature and oxygen transfer rate) permit such growth. The term "basal medium" refers to a medium which promotes the growth of many types of microorganisms which do not require any special nutrient supplements. Most basal media generally comprise of four basic chemical groups: amino acids, carbohydrates, inorganic salts, and vitamins. A basal medium generally serves as the basis for a more complex medium, to which supplements such as serum, buffers, growth factors, lipids, and the like are added. In one aspect, the growth medium may be a complex medium with the necessary growth factors to support the growth and expansion of the cells of the disclosure while maintaining their self-renewal capability. Examples of basal media include, but are not limited to, Eagles Basal Medium, Minimum Essential Medium, Dulbecco's Modified Eagle's Medium, Medium 199, Nutrient Mixtures Ham's F-10 and Ham's F-12, McCoy's 5A, Dulbecco's MEM/F-12, RPMI 1640, and Iscove's Modified Dulbecco's Medium (IMDM).

Culture Period

[0167] In the context of this disclosure, "culturing" cells for specified periods of time (e.g. at least 10 weeks) refers to a time period wherein day zero or "day 0" is the time point at which the cells are transferred to the culture vessel. The culture vessel may be a flask, for example the standard T-175 cell culture flask. Typically, the culture vessel is a multi-compartment bioreactor such as the Integra CELLline bioreactor, and day zero is the day on which the stem cells are transferred into the bioreactor. Accordingly, cells "that have been cultured for at least 10 weeks" refers to cells that have been cultured for at least 10 weeks following transfer into the culture vessel. In this 10 week period, the cells are not passaged or subcultured, i.e. they are not transferred to a new culture vessel. Optionally, cells can be removed from the culture vessel during the culture period, typically for sampling, but this does not change the cells that remain in the culture vessel, which have been in that culture vessel since day 0.

[0168] In one embodiment, as described in Example 10, on day zero approximately 15×10^6 CTX0E03 cells in a total of 15ml of complete growth medium are introduced into the cell compartment of the CeLLline bioreactor, followed by the addition of a further 460ml of complete growth medium to the cell compartment.

Pharmaceutical Compositions

[0169] The neural stem cell exosome of the invention is useful in therapy and can therefore be formulated as a pharmaceutical composition. A pharmaceutically acceptable composition typically includes at least one pharmaceutically acceptable carrier, diluent, vehicle and/or excipient in addition to the exosomes of the invention. An example of a suitable carrier is Ringer's Lactate solution. A thorough discussion of such components is provided in Gennaro (2000) Remington: The Science and Practice of Pharmacy. 20th edition, ISBN: 0683306472.

[0170] The phrase "pharmaceutically acceptable" is employed herein to refer to those compounds, materials, compositions, and/or dosage forms which are, within the scope of sound medical judgment, suitable for use in contact with the tissues of human beings and animals without excessive toxicity, irritation, allergic response, or other problem or complication, commensurate with a reasonable benefit/risk ratio.

[0171] The composition, if desired, can also contain minor amounts of pH buffering agents. The carrier may comprise storage media such as Hypothermosol®, commercially available from BioLife Solutions Inc., USA. Examples of suitable pharmaceutical carriers are described in "Remington's Pharmaceutical Sciences" by E W Martin. Such compositions will contain a prophylactically or therapeutically effective amount of a prophylactic or therapeutic microparticle preferably in purified form, together with a suitable amount of carrier so as to provide the form for proper administration to the subject. The formulation should suit the mode of administration. In a preferred embodiment, the pharmaceutical compositions are sterile and in suitable form for administration to a subject, preferably an animal subject, more preferably a mammalian subject, and most preferably a human subject.

[0172] The pharmaceutical composition may be in a variety of forms. These include, for example, semisolid, and liquid dosage forms, such as lyophilized preparations, liquid solutions or suspensions, injectable and infusible solutions. The

pharmaceutical composition is preferably injectable. A particular advantage of the exosomes of the invention is their improved robustness compared to the stem cells from which they are obtained; the exosomes can therefore be subjected to formulation, such as lyophilisation, that would not be suitable for stem cells.

[0173] Pharmaceutical compositions will generally be in aqueous form. Compositions may include a preservative and/or an antioxidant.

[0174] To control tonicity, the pharmaceutical composition can comprise a physiological salt, such as a sodium salt. Sodium chloride (NaCl) is preferred, which may be present at between 1 and 20 mg/ml. Other salts that may be present include potassium chloride, potassium dihydrogen phosphate, disodium phosphate dehydrate, magnesium chloride and calcium chloride.

[0175] Compositions may include one or more buffers. Typical buffers include: a phosphate buffer; a Tris buffer; a borate buffer; a succinate buffer; a histidine buffer; or a citrate buffer. Buffers will typically be included at a concentration in the 5-20mM range. The pH of a composition will generally be between 5 and 8, and more typically between 6 and 8 e.g. between 6.5 and 7.5, or between 7.0 and 7.8.

[0176] The composition is preferably sterile. The composition is preferably gluten free. The composition is preferably non-pyrogenic.

[0177] In a typical embodiment, the microparticles are suspended in a composition comprising 6-hydroxy-2,5,7,8-tetramethylchroman-2-carboxylic acid (Trolox®), Na⁺, K⁺, Ca²⁺, Mg²⁺, Cl⁻, H₂PO₄⁻, HEPES, lactobionate, sucrose, mannitol, glucose, dextran-40, adenosine and glutathione. Typically, the composition will not include a dipolar aprotic solvent, e.g. DMSO. Suitable compositions are available commercially, e.g. HypoThermasol®-FRS. Such compositions are advantageous as they allow the microparticles to be stored at 4°C to 25°C for extended periods (hours to days) or preserved at cryothermic temperatures, i.e. temperatures below -20°C. The microparticles may then be administered in this composition after thawing.

[0178] The pharmaceutical composition can be administered by any appropriate route, which will be apparent to the skilled person depending on the disease or condition to be treated. Typical routes of administration include intravenous, intra-arterial, intramuscular, subcutaneous, intracranial, intranasal or intraperitoneal. For treatment of a disorder of the brain, one option is to administer the microparticles or miRNA intra-cerebrally, typically to the site of damage or disease.

[0179] The microparticles or miRNA will be administered at a therapeutically or prophylactically-effective dose, which will be apparent to the skilled person. Due to the low or non-existent immunogenicity of the microparticles, it is possible to administer repeat doses without inducing a deleterious immune response.

Therapeutic uses

[0180] The microparticles and miRNA of the disclosure are useful in the treatment or prophylaxis of disease. Accordingly, the disclosure includes a method of treating or preventing a disease or disorder in a patient using a microparticle or miRNA of the disclosure. The term "patient" includes human and other mammalian subjects that receive either prophylactic or therapeutic treatment.

[0181] The Examples below demonstrate that neural stem cell exosomes have been identified that *inhibit* cell migration. Inhibition of migration has been observed in human fibroblasts. This inhibition is particularly surprising because, as also shown in the examples and as described in PCT/GB2013/050879, neural stem cell microparticles have previously been shown to *stimulate* fibroblast migration. Inhibition of migration has also been observed in glioblastoma cells.

[0182] The exosomes of the invention are used in the therapy of glioma, for example glioblastoma. In a further embodiment, the exosomes of the invention treat the cancer by inhibiting migration of the cancer cells. In yet a further embodiment, the microparticles of the invention treat the cancer by inducing differentiation of cancer cells, typically differentiation of a nestin-positive cancer cell. In another aspect described herein, the microparticles may treat the cancer by inducing or enhancing an immune response against the cancer cells. Glioma is a CNS cancer, so the immune response typically comprises the activation and/or proliferation of glial cells such as microglia.

[0183] The cancer may be a solid tumour cancer, for example a sarcoma or carcinoma. The solid tumour cancer may also be a solid lymphoma. Exemplary solid tumour cancers include breast cancer, lung cancer, prostate cancer, bowel cancer, renal cancer, hepatic cancer, pancreatic cancer, cervical cancer, testicular cancer, gastric (stomach) cancer, uterine cancer, ovarian cancer, cancers of the head and neck, mouth cancer, thyroid cancer, oesophagus cancer, brain cancer including glioma (e.g. glioblastoma) and meningioma, Kaposi's sarcoma, Castleman's disease, cutaneous T-cell lymphoma (CTCL), cutaneous B-cell lymphoma, and skin cancer such as basal cell carcinoma, squamous cell carcinoma and melanoma.

[0184] In one embodiment, the solid tumour cancer is breast cancer, typically ductal carcinoma in situ, lobular carcinoma in situ, invasive ductal carcinoma, invasive lobular carcinoma, inflammatory breast cancer or Paget's disease.

[0185] In another embodiment, the solid tumour cancer is lung cancer, typically squamous cell carcinoma, adenocarcinoma or large cell carcinoma, or a small cell lung cancer.

[0186] In a further embodiment, the solid tumour cancer is prostate cancer, typically prostate adenocarcinoma.

[0187] In a further embodiment, the solid tumour cancer is skin cancer, typically a basal cell carcinoma, squamous cell carcinoma or melanoma.

[0188] The cancer may be a liquid tumour, which is typically a tumour of the blood, bone marrow, or lymph nodes. Such cancers include leukemia, lymphoma and myeloma. Exemplary liquid tumours include acute lymphoblastic leukemia, acute myelogenous leukemia (AML), multiple myeloma, Hodgkin's lymphoma and non-Hodgkins lymphoma.

[0189] The cancer is a glioma. An exemplary glioma is a glioblastoma, which may be a giant cell glioblastoma or a gliosarcoma. The *in vivo* xenograft pilot data in the Examples demonstrate trends in the treatment of glioblastoma.

[0190] The Examples below demonstrate that exosomes of the invention (in the case of Example 4, exosomes isolated from proliferating CTX0E03 cell culture) reduce the expression of nestin on tumour cells. Accordingly, in one embodiment, the cancer to be treated is nestin-positive. Nestin-positive cancers include melanoma, breast cancer, CNS cancers such as glioma and typically glioblastoma, pancreatic cancer, gastrointestinal stromal tumours (GISTs), dermatofibrosarcoma protuberances, thyroid tumours and prostate cancer (see, for example, Ishiwata et al World J Gastroenterol. 2011 January 28; 17(4):409-418). The nestin-positive breast cancer is typically "triple negative, nestin positive" breast cancer (ER α /PR $^-$ /Her2 $^-$ /Nestin $^+$). Triple negative breast cancer is an aggressive disease, recurring and metastasizing more often than other kinds of breast cancer, and treatments for this are urgently needed. The effectiveness of microparticles of the invention in treating this cancer can readily be tested *in vivo* using a triple negative breast cancer mouse model, for example as described by Kaur et al, BMC cancer 2012m 12:120. *In vivo* models for other cancers exist and can be used to test the effectiveness of microparticles of the invention; for example, xenograft models of melanoma (e.g. Rofstad Br. J. Cancer (1994), 70, 804-812) and glioblastoma (e.g. Jacobs et al, ASN Neuro. 2011; 3(3); 2011).

[0191] Nestin is also reported to be expressed in endothelial cells involved in angiogenesis (Mokry et al, Stem Cells Dev. 2004; 13:658-664) and so the ability of microparticles of the invention to reduce nestin expression provides a further mechanism to inhibit angiogenesis.

[0192] Microparticles of the disclosure may also be used to treat or prevent metastatic cancers, for example metastasis of each of the cancers listed above.

[0193] Fibroblasts are known to play a role in angiogenesis during tumour formation. Without being bound by theory, it is thought that this is mediated in part by a paracrine mechanism wherein factors secreted by the fibroblasts, including Fibroblast Growth Factor (FGF), act on endothelial cells in the nascent or growing blood vessel. Therefore, inhibiting the migration of fibroblasts is expected to inhibit angiogenesis. Accordingly, the microparticles of the disclosure may be used as an anti-angiogenic therapy, i.e. in the therapy of unwanted, deleterious or undesirable angiogenesis. In one embodiment, the unwanted or undesirable angiogenesis is a component or a precursor of a solid tumour, typically a cancerous solid tumour. In this embodiment, the microparticles are used in the therapy of the tumour by preventing, inhibiting or reducing angiogenesis in the tumour. Typically, the solid tumour that is treated by targeting angiogenesis is one of the tumours described above, for example a sarcoma or carcinoma. The solid tumour cancer in this embodiment may also be a solid lymphoma. Exemplary solid tumour cancers that can be treated by targeting the angiogenic component of the tumour include breast cancer, lung cancer, prostate cancer, bowel cancer, renal cancer, hepatic cancer, pancreatic cancer, cervical cancer, testicular cancer, gastric (stomach) cancer, uterine cancer, ovarian cancer, cancers of the head and neck, mouth cancer, thyroid cancer, oesophagus cancer, brain cancer including glioma (e.g. glioblastoma) and meningioma, Kaposi's sarcoma, Castleman's disease, cutaneous T-cell lymphoma (CTCL), cutaneous B-cell lymphoma, and skin cancer such as basal cell carcinoma, squamous cell carcinoma and melanoma.

[0194] In one embodiment, the microparticles and compositions containing them are not used for immune modulation. In one embodiment, the therapy is not related to immunomodulation.

[0195] The disclosure also provides a method for treating or preventing a disease or condition comprising administering an effective amount of the microparticle of the disclosure, thereby treating or preventing the disease. Typically, the disease or condition is as identified above.

[0196] In one embodiment, the microparticles for use in therapy are isolated from NSCs (typically CTX0E03 cells) that have been cultured (typically in a multi-compartment bioreactor) for at least 10 weeks, typically at least 11 weeks, at least 12 weeks, at least 13 weeks, at least 14 weeks or at least 15 weeks. Optionally, the NSCs have been cultured for no more than 20 weeks, e.g. between 10 and 20 weeks, between 11 and 20 weeks, between 12 and 20 weeks, between 13 and 20 weeks, between 14 and 20 weeks or between 15 and 20 weeks. Typically, the microparticles are exosomes. In the examples, microparticles produced according to this embodiment are shown to inhibit fibroblast migration and induce or enhance tumour destruction by the immune system.

[0197] The observed increased efficacy of exosomes isolated from NSCs (CTX0E03 cells) that have been cultured (in a multi-compartment bioreactor) for 6 weeks correlates with the observed reduction in size of the exosomes to around 70nm diameter, which also occurred after culturing the cells for 6 weeks. Accordingly, in one embodiment exosomes isolated from NSCs (typically CTX0E03 cells) having a diameter less than 100nm, typically less than 80nm, for example around 70nm diameter, are used in therapy as described above.

[0198] In another embodiment, the microparticles for use in therapy are isolated from proliferating NSCs (typically CTX0E03 cells) that have been cultured in a standard culture vessel such as a T-175 flask, or have been cultured in a

multi-compartment bioreactor for 4 weeks or less, 3 weeks or less, 2 weeks or less, or 1 week or less e.g. exosomes isolated on day 0 of the multi-compartment culture. These cells are typically passaged when sub-confluent, are positive for a stem cell marker (e.g. nestin) and negative for markers of differentiated cells (e.g. GFAP or DCX). These exosomes may have a diameter greater than 100nm. In the examples, microparticles produced according to this embodiment are shown to inhibit cancer cell migration and induce tumour cell differentiation.

[0199] In prophylactic applications, pharmaceutical compositions or medicaments are administered to a patient susceptible to, or otherwise at risk of, a particular disease in an amount sufficient to eliminate or reduce the risk or delay the outset of the disease. In therapeutic applications, compositions or medicaments are administered to a patient suspected of, or already suffering from such a disease in an amount sufficient to cure, or at least partially arrest, the symptoms of the disease and its complications. An amount adequate to accomplish this is defined as a therapeutically-or pharmaceutically-effective dose. In both prophylactic and therapeutic regimes, agents are typically administered in several dosages until a sufficient response has been achieved. Typically, the response is monitored and repeated dosages are given if the response starts to fade.

[0200] The microparticles of the disclosure may optionally be combined with a stem cell to provide a combination therapy. The stem cell is optionally the stem cell from which the microparticle is derived, e.g. if the microparticle is an exosome from a CTX0E03 cell, then the stem cell for use in combination therapy may be a CTX0E03 cell, typically but not necessarily cultured for the same period of time as the cells from which the microparticles were derived. A stem cell and microparticle can optionally be (i) administered together in a single pharmaceutical composition, (ii) administered contemporaneously or simultaneously but separately, or (iii) administered separately and sequentially, e.g. stem cell followed by microparticle, or microparticle followed by stem cell. When the stem cell and microparticle are administered separately and sequentially, the duration between the administration of the cell and microparticle may be one hour, one day, one week, two weeks or more.

[0201] In one embodiment, a prophylactic therapy induces tolerance, typically immunotolerance, in a host that is to receive the stem cells from which the microparticle is derived. In one embodiment, the administration of one or more doses of microparticles of the disclosure to a patient, prior to administration of a stem cell therapy, can be used to reduce the risk of an adverse immune response, i.e. "rejection", of the stem cell therapy. In another embodiment, tolerance to the stem cells can be increased by administering stem cells together with microparticles of the disclosure, as discussed above.

Effective doses of the compositions of the present invention, for the treatment of the above described conditions vary depending upon many different factors, including means of administration, target site, physiological state of the patient, whether the patient is human or an animal, other medications administered, and whether treatment is prophylactic or therapeutic. Usually, the patient is a human.

[0202] The CTX0E03 cell line has been shown to be effective in treating stroke, peripheral arterial disease, brain damage such as motor, sensory and/or cognitive deficit, and psychiatric disorders. The cells are currently being tested in a clinical trial for treatment of disabled stroke patients (Clinicaltrials.gov Identifier: NCT01151124). WO-A-2012/004611 describes the use of the CTX0E03 cells in treating psychiatric disorders including unipolar and bipolar depression, schizophrenia, obsessive compulsive disorder, autism and autistic syndrome disorders.

[0203] As used herein, the terms "treat", "treatment", "treating" and "therapy" when used directly in reference to a patient or subject shall be taken to mean the amelioration of one or more symptoms associated with a disorder, or the prevention or prophylaxis of a disorder or one or more symptoms associated with a disorder. The disorders to be treated include, but are not limited to, cancer, fibrosis, rheumatoid arthritis, atherosclerosis, and other diseases involving deleterious cell migration. Amelioration or prevention of symptoms results from the administration of the microparticles of the disclosure, or of a pharmaceutical composition comprising these microparticles, to a subject in need of said treatment.

Tracing administered cells and microparticles *in vivo*

[0204] The present disclosure provides a distinct marker profile for microparticles produced by neural stem cells. It is therefore possible to detect the presence of these microparticles *in vivo*, by testing a sample obtained from a patient and determining whether the marker profile in the sample matches that of the microparticles. If the sample profile matches the profile of the microparticles described herein, then this confirms the presence of the microparticles. This can be used to detect not only the presence and/or biodistribution of the microparticles themselves, but also the presence of stem cells producing the microparticles. This is particularly useful when detecting whether a stem cell administered *in vivo* has engrafted into the host tissue, and/or has migrated, for example in ADME(T) studies.

[0205] Detection of the microparticles *in vivo* can be used to monitor the course of a treatment wherein microparticles or stem cells are administered to a patient. Determining the presence, absence or amount of microparticles or cells producing microparticles of the disclosure in a patient allows the dosage regime to be altered accordingly, e.g. to increase or decrease the dose as required to provide an effective amount of microparticles or stem cells *in vivo*.

Methods of producing microparticles

[0206] Microparticles are isolated from stem cell conditioned media. The "conditioned medium" (CM) may be a growth medium for stem cells, which has been used to culture a mass culture of stem cells for at least about 12 hours, at least about 24 hours, at least about 48 hours or least about 72 hours, typically up to 168 hours (7 days), removed and sterilized by any suitable means, preferably by filtration, prior to use, if required.

[0207] Microparticles that are able to inhibit fibroblast cell migration have been isolated from stem cells that have been cultured for at least 10 weeks. Accordingly, one way to produce microparticles that are able to inhibit cell migration is to culture the cells in a multi-compartment bioreactor for at least about 10 weeks before the microparticles are harvested, typically at least 11 weeks, at least 12 weeks, at least 13 weeks, at least 14 weeks, at least 15 weeks, and optionally no longer than 20 weeks. Example 10 describes a typical culture protocol using a CeLLine bioreactor.

[0208] Microparticles that are able to inhibit glioblastoma cell migration have been isolated from proliferating stem cells that have been cultured for 4 weeks or less. Accordingly, one way to produce microparticles that are able to inhibit cell migration is to culture the cells so that they are able to proliferate, for example by culturing in a T-175 flask, or in a multi-compartment bioreactor for 4 weeks or less, 3 weeks or less, 2 weeks or less, or 1 week or less e.g. exosomes isolated on day 0 of the multi-compartment culture.

[0209] Typically, microparticles may be harvested from a multi-compartment, e.g. two-compartment, bioreactor which allows the cell culture, and hence the conditioned media, to be maintained for longer periods of time, for example more than 10 weeks, at least 11 weeks, at least 12 weeks, at least 13 weeks, at least 14 weeks, at least 15 weeks, and optionally no longer than 20 weeks.. The system maintains the cells and secreted microparticles within a small cell compartment (approximately 15ml) which is separated from a larger reservoir of medium by a 10kDa semi-permeable membrane. This allows the efficient removal of metabolic waste products while effectively maintaining an extremely high cell density to maximize microparticle production. Example 14, and Figures 9 and 10, demonstrate that use of a two-compartment bioreactor results in a much higher yield of microparticles than is obtained when a standard cell culture flask (T175 flask) is used.

[0210] The microparticles may be separated from other media components based on molecular weight, size, shape, hydrodynamic radius, composition, charge, substrate-ligand interaction, absorbance or scattering of electromagnetic waves, or biological activity. In one embodiment, the conditioned media is filtered using a filter of appropriate size to separate the desired microparticle, for example a 100K MWCO filter. Optionally, the stem cell-conditioned medium is concentrated prior to the isolation of the microparticles by subjecting the concentrated NSC-conditioned medium to size exclusion chromatography. The UV absorbant fractions can then be selected for isolation of the microparticles of interest.

[0211] Different microparticles can be isolated from the media by using different isolation techniques and parameters. For example, exosomes have a vesicle density of 1.13-1.19 g/mL and can be isolated by differential centrifugation and sucrose gradient ultracentrifugation at 100,000-200,000g. Microvesicles can be isolated by filtration (100K MWCO) and differential centrifugation at 18,000-20,000g. Membrane particles have a density of 1.04-01.07 g/ml and Exosome-like vesicles have a density of 1.1 g/ml.

[0212] A typical production method comprises: culturing stem cells to produce conditioned media; removing cell debris by centrifugation at 1500 rpm; isolating microvesicles (<1000kDa) by ultrafiltration through a 100K MWCO filter or isolating exosomes (30-100nm) by ultracentrifugation at 120,000g; followed by quantification using a BCA protein assay.

Conditionally immortalised stem cells as producer cells for microparticles

[0213] In one aspect of the disclosure, conditionally immortalised stem cells are used to produce microparticles such as microvesicles and/or exosomes. These conditionally immortalised stem cells are typically neural stem cells, but may be a stem cell of any type, for example a haematopoietic stem cell or a mesenchymal stem cell. A method of producing stem cell microparticles is therefore provided, comprising the steps of culturing conditionally-immortalised stem cells and harvesting the microparticles that are produced by the cells, as described above. Conditional immortalisation of stem cells is known in the art, as described above. For the avoidance of doubt, this method is not limited to the use of neural stem cells.

[0214] When the stem cell used to produce microparticles is a neural stem cell, it may be any of the neural stem cells described herein, for example the CTX0E03 conditionally-immortalised cell line which is clonal, standardised, shows clear safety *in vitro* and *in vivo* and can be manufactured to scale thereby providing a unique resource for stable exosome production. Alternatively, the neural stem cells may be neural retinal stem cell lines, optionally as described in US 7514259.

[0215] When the stem cell used to produce microparticles is a mesenchymal stem cell, it may optionally be a conditionally-immortalised adipose-derived stem cell ("ADSC") or a conditionally-immortalised version of the mesenchymal stem cells described in WO-A-2009/105044; these cells are CD29+, CD44+, CD49a+/e+, CD105+, CD166+, CD34-, CD45-.

Methods of inducing microparticle secretion

[0216] The inventors have found that it is possible to increase the production of microparticles by stem cells. This finding, which is not limited to neural stem cells and can be used for the production of microparticles from any stem cell, allows for an improved yield of microparticles to be obtained from a stem cell culture.

[0217] A first technique to increase the production of microparticles by the stem cells is to treat the stem cells with one or more of TGF- β , IFN- γ or TNF- α , typically at between 1 and 25ng/ml e.g. 10ng/ml, for between 12 to 96 hours prior to the removal of conditioned media.

[0218] As explained in Example 8 below, the frequency of the occurrence of multivesicular bodies (MVBs) was observed to be altered by the presence of TGF- β , IFN- γ or TNF- α (10ng/ml). The frequency was highest in the presence of TGF- β , followed by IFN- γ , followed by TNF- α . Therefore, adding one or more of TGF- β , IFN- γ or TNF- α to the stem cell culture medium will stimulate the production of microparticles by the cells. The microparticles can then be harvested, by separating the microparticles from other components as described above.

[0219] A second technique to increase the production of microparticles by the stem cells is to culture the cells under hypoxic conditions. Culturing cells under hypoxic conditions is well-known to the skilled person, and involves culturing the cells in an atmosphere that has less than atmospheric level of O₂, i.e. less than 21% O₂. This is typically achieved by placing the cells in an incubator that allows oxygen levels to be changed. Hypoxic culture typically involves culturing in an atmosphere containing less than 10% O₂, more typically 5% or less O₂, for example 4% or less, 3% or less, 2% or less, or 1% or less O₂.

[0220] The inventors have also realised that co-culturing a stem cell with a different cell type can alter the production of microparticles by the stem cell. The different cell type may be a non-stem cell, i.e. a terminally differentiated cell type. Typically, the different cell type is one with which the stem cell would interact *in vivo*. In one embodiment, neural stem cells are co-cultured with epithelial cells such as endothelial cells, typically Human Umbilical Vein Endothelial Cells (HUVEC). It has been observed that *in vivo*, NSCs and the vasculature interact, with proliferating NSCs being localized in close proximity or adjacent to blood vessels. Receptor tyrosine kinase activation and signal protein secretion has also been observed to be upregulated when NSCs are co-cultured with endothelial cells, again indicating that the vasculature modulates the proliferation capacity of NSCs.

[0221] Therefore, culturing a stem cell with a different cell type may improve the amount of microparticles produced and/or may refine the content of the microparticles, typically so that the microparticles produced by the stem cells are further biased towards a state of inhibition of cell migration. Accordingly, microparticles produced by stem cells that have been co-cultured with other cells, e.g. NSCs co-cultured with endothelial cells, are advantageous. These microparticles may be obtained by isolation from the co-cultured stem-cell conditioned media, as described herein.

[0222] Surprisingly, the present inventors have realised that the amount of microparticles produced by stem cells can be increased greatly simply by culturing stem cells in a multi-compartment bioreactor. This finding is not limited to neural stem cells and applies generally to the culture of all stem cells. Accordingly, one aspect of the disclosure provides a method of producing microparticles from stem cells that have been cultured in a multi-compartment bioreactor. The cells from which the microparticles are harvested have typically been cultured for at least one week, typically at least 8, 9, 10, 11, 12, 13 or 14 days, for example 15 days, 16 days, 17 days, 18 days, 19 days, 20 days, 21 days or more, for example at least three weeks, four weeks, five weeks, six weeks or more. To produce microparticles that inhibit cell migration, the cells from which the microparticles are harvested have typically been cultured for more than ten weeks. It can be seen from Figure 10 that the increase in microparticle production, week on week, is not merely additive but is exponential. The prolonged culture typically has been observed in the Integra Celline system two-compartment bioreactor (commercially available from Integra Biosciences AG, Zizers, Switzerland) but the findings are not limited to this specific multi-compartment bioreactor; any multi-compartment bioreactor can be used. This culture method can be used to produce microparticles from any stem cell type, including but not limited to neural stem cells and mesenchymal stem cells.

Method for screening total RNA composition of conditioned medium

[0223] Following centrifugation (5 min at 1500 rpm), microparticles are collected from conditioned medium through filtration (0.02-0.2 μ m, or 100K MWCO). Total RNA is obtained using trizol based extraction followed by purification using Qiagen RNeasy mini kit. The extract in water has a 260:280 nm absorbance suggesting that it may be RNA. Total RNA is retro-transcribed with either a protocol suitable for mRNA (Superscript II RT, Invitrogen) or miRNA (mScript RT kit, Qiagen). Validation of mRNA and miRNA presence is proven by qRT-PCR using primers for ATP5B and YWHAZ for mRNA, and U6B and 15a for miRNA housekeeping genes respectively. The RNA may further be assessed by a generic gene expression analysis assay such as an array (micro array or PCR based array), and sequencing.

[0224] The invention is further described with reference to the following non-limiting examples.

Examples

Example 1: NSC exosomes that inhibit cell migration.

[0225] A transwell assay was used to study the migratory response of human dermal fibroblasts to different populations of exosomes. Experiments were performed in triplicate. 200,000 human dermal fibroblast cells ("FBs") were placed on the upper layer of a cell permeable membrane (8 μ m pore size; 24-well plate) and a solution (basal medium) containing or lacking 20 μ g/ml exosomes was placed in contact with the underside of the cell permeable membrane (Figure 1, top panel). The exosomes were collected from CTX0E03 cells cultured for 0 weeks ("0") or 11 weeks ("11") in an Integra

CellLine AD1000 multi-chamber bioreactor. Following an incubation period (6 or 24 hours; control: 0 hours), the human dermal fibroblast cells that migrated through the membrane were stained (using a fluorescent-dye conjugated anti-actin antibody and Hoechst Fluorescent Stain for nuclei) and counted (six random microscope fields per sample) as an indicator of the cells' migratory response to exosomes.

[0226] Figure 1 (lower panel) and Figure 2 show that exosomes isolated from a proliferating CTX0E03 culture ("0") significantly promote migration of human dermal fibroblasts compared to medium lacking exosomes ("basal"), both after a 6 hour and after a 24 hour incubation period. In contrast, exosomes isolated from a more differentiated CTX0E03 culture ("11") significantly abrogate migration of human dermal fibroblasts compared to medium lacking exosomes ("basal"), both after a 6 hour and after a 24 hour incubation period.

[0227] It can be seen that cell migration is increased in the presence of exosomes from 0-week NSCs but decreased in the presence of exosomes from 11-week NCSs, compared to control ("basal").

[0228] In summary, NSC microparticles have been identified that significantly abrogate cell migration.

[0229] These data show that neural stem cell microparticles can stimulate or inhibit cell migration. This is surprising and useful in applications where either stimulating (e.g. wound healing) or inhibiting (e.g. cancer, fibrosis, rheumatoid arthritis, atherosclerosis) cell migration is desired. The involvement of fibroblasts in angiogenesis also makes the microparticles that inhibit fibroblast migration useful in applications where inhibition of angiogenesis is desired. Angiogenesis is involved in tumour formation, survival and metastasis. These data therefore demonstrate potential for the exosomes of the invention to treat many types of cancer.

Example 2: Exosomes isolated from the medium of NSCs cultured for 2 or 6 weeks promote fibroblast migration.

Method - Wound closure/ scratch assay

[0230]

- Seed 0.25×10^6 NHDF (normal human dermal fibroblasts) per well of a 12 well plate and allow to become confluent (24 hours)
- Remove growth factors for 24hrs
- Remove cells (scratch) and incubate with exosomes/conditioned media
- Image effected area over 48hrs
- Estimate area using Image J

Results

[0231]

Table 2 - Wound closure/scratch assay representing the migration activity of normal human dermal fibroblasts (NHDF) cultured in CTX0E03 conditioned media or upon the addition of purified exosomes.

	Wound closure (%)		
	0h	24h	48h
CTX0E03 conditioned media	0%	100%	
2ug/ml exosomes	0%	95.4%	100%
Control	0%	48.1%	49.7%

[0232] Wound closure was calculated as the area covered by cells in relation to the initial wound area, as determined at 0h. Wound closure is expressed as the percentage of the initial wound area at time 0h. These data are also shown,

photographically, in Figure 3A. The figure shows that, in contrast to exosomes from 11-week NSCs as described in Figures 1 and 2, exosomes from 2-week NSCs stimulate cell migration.

[0233] Figure 3B shows that 10 μ g CTX0E03 exosomes significantly increase wound closure (as determined in the HDNF scratch/migration assay) after 72 hours, compared to basal conditions (without exosomes).

[0234] Further experiments confirmed that exosomes purified (by ultracentrifugation; quantified by BCA protein assay; characterised as >99% positive for CD63 and CD81 and having a greater expression level of Alix compared to the corresponding microparticle fraction) from all time points (weeks 2-6) during continuous culture (using Integra CELLline bioreactors in the presence of growth factors and 4OHT) significantly enhanced fibroblast migration and wound healing, with a peak response between 5-10 μ g/ml compared to basal conditions. Figure 3C shows the % healed areas for basal conditions, 2 μ g/ml exosomes, 6 μ g/ml exosomes, 20 μ g/ml exosomes and an LSGS (low serum growth supplement) positive control. The top panel of Figure 3C shows exosomes isolated from CTX0E03 cells cultured for 2 weeks in the Integra Cellline system and the bottom panel of Figure 3C shows exosomes isolated from CTX0E03 cells cultured for 6 weeks in the Integra Cellline system. These data show that all doses of all tested NSC exosomes provide increased healing compared to basal conditions, with % healing approaching the positive control (LSGS) after 72 hours.

[0235] The data in Figure 3C also show that the exosomes isolated from NSCs cultured for 6 weeks cause faster healing (than 2 week exosomes), with the % healed approaching 100% after only 48 hours, for all doses.

[0236] Figure 3D shows the results of an *in vivo* injection wound assay in a mouse, confirming that CTX0E03 cells stimulated wound healing to a statistically-significant degree *in vivo*. This is a simple *in vivo* bioassay which can be used to confirm the efficacy of microparticles *in vivo*.

[0237] **Conclusion** Exosomes released from the human neural stem cell line CTX0E03 enhance fibroblast migration in an *in vitro* model of wound healing, suggesting that exosomes may contribute to the mechanisms by which hNSCs promote repair. Exosomes isolated from cells cultured for 6 weeks show improved wound healing efficacy *in vitro*, compared to exosomes isolated from cells cultured for 2 weeks.

Example 3: Glioblastoma engraftment assay - Destruction of tumour cells

[0238] U373 glioblastoma cells were pre-treated *in vitro* for 24 hours with exosomes isolated from CTX0E03 cells cultured for 11 weeks in an Integra CeLLine bioreactor before implantation into the striatum of Balb-C mice brains.

[0239] As shown in Figure 19, the exosome-treated glioblastoma cells did not engraft into the striatum. Histopathology demonstrated the presence of necrotic U373 cell bodies at the site of implantation and evidence of gliosis - a host cellular immune response.

[0240] These data suggest utility of these exosomes in the treatment of cancer, by promoting the destruction of a tumour by the immune system, particularly a tumour of the CNS such as a glioblastoma.

Example 4: Glioblastoma engraftment assay - Differentiation of tumour cells

[0241] U373 glioblastoma cells were pre-treated *in vitro* for 24 hours with exosomes isolated from standard CTX0E03 cell culture ("exosome 0" - proliferating cells, cultured in an Integra CeLLine bioreactor for less than 24 hours) before implantation into the striatum of Balb-C mice brains. Marker expression was then observed after 24 hours.

[0242] As shown in Figure 20, the exosome-treated glioblastoma cells demonstrated a reduction in nestin expression 24 hours post implantation into the striatum of Balb-C mice. Nestin is a stem cell marker; cancer stem cells drive tumourigenesis, are linked with metastasis, high grade and poor prognosis. The treatment of cancer by inducing cellular differentiation is particularly attractive because the therapy can be target-cell specific (i.e. will only target the undifferentiated, malignant, cells) and likely less toxic than standard chemotherapies.

[0243] These data suggest utility of these exosomes in the treatment of cancer, by inducing differentiation of the cancer cells, typically for treating a nestin-positive cancer and particularly a tumour of the CNS such as a glioblastoma.

Example 5: *In vitro* glioblastoma differentiation assay - Differentiation of tumour cells

[0244] U373 glioblastoma cells were cultured for 24 hours in the presence of: (i) basal medium; (ii) +20 μ g exosomes isolated from standard CTX0E03 cell culture ("exosome 0" - proliferating cells, cultured in an Integra CeLLine bioreactor for less than 24 hours); or (iii) +20 μ g exosomes isolated from CTX0E03 cells cultured for 11 weeks in an Integra CeLLine bioreactor ("exosome 11"). The U373 cells were then stained for the presence of Nestin (a stem cell marker) and GFAP (an astrocyte marker of a differentiated cell).

[0245] As shown in Figure 21, exosome 0 promoted differentiation of the glioblastoma cells *in vitro*. The exosome 0 treated cells appeared morphologically differentiated, with the presence of long processes. Additionally, these cells expressed Glial fibrillary acidic protein (GFAP), a marker of differentiated astroglial cells. These *in vitro* data agree with the *in vivo* data immediately above. As noted above, more differentiated (less malignant) glioblastoma tumours are linked

with more favourable prognosis. These data further suggest utility of these exosomes in the treatment of cancer, by inducing differentiation of the cancer cells, particularly a tumour of the CNS such as a glioblastoma.

[0246] In contrast, the exosomes isolated from CTX0E03 cells cultured for 11 weeks ("exosome 11") promoted "stemness" in the glioblastoma cells *in vitro*, demonstrated by nestin expression and proliferation. However, in the *in vivo* assay above, these exosomes were observed to promote destruction of the tumour cells.

Example 6: Glioblastoma migration assays

[0247] Three separate *in vitro* transmembrane migration assays have demonstrated that treatment of glioblastoma (U373) cells with exosomes isolated from standard CTX0E03 cell culture ("exosome 0" - proliferating cells, cultured in an Integra CeLLine bioreactor for less than 24 hours) significantly reduces their migration towards a positive chemoattractant (Foetal Bovine Serum). These assay results are shown in Figure 22.

[0248] Glioblastoma cells were seeded on one side of a porous filter membrane. These cells were either seeded together with 20 μ g/ml CTX0E03 "exosome 0" (Figures 22A and 22C) or have been pre-treated with 10 μ g/ml CTX0E03 "exosome 0" for 24 hours (Figure 22B). Medium containing a 10%FBS was placed on the opposing (lower) side. After a 24 hour incubation period, the membrane was fixed and stained to reveal migrated cells (e.g. cell nuclei), which were counted microscopically.

[0249] These data show that exosomes of the invention (isolated from standard "week 0" CTX0E03 cells) are able to reduce migration of glioblastoma cells. Glioblastomas are the most common and malignant brain tumors of the central nervous system and exhibit high invasive capacity, which hinders effective therapy. Therefore, therapeutics that can inhibit glioma cell migration and invasion are highly desirable. These data demonstrate the utility of neural stem cell exosomes in the treatment of cancer, typically a glioblastoma, by reducing tumour migration/invasion.

Summary: Treatment of cancer using neural stem cell exosomes

[0250] The data provided above indicate therapeutic utility in the treatment of cancer using exosomes produced by neural stem cells, by one or more of: reducing tumour migration/invasion (exosome 0, glioblastoma assay); inducing tumour differentiation (exosome 0, glioblastoma assay); promoting tumour destruction (exosome 11, glioblastoma transplant); or inhibiting angiogenesis (exosome 11, fibroblast assay).

Example 7: Preparation of neural stem cells and neural stem cell microparticles for visualisation by electron microscopy.

Method

Embedding CTX0E03 cells for electron microscopy

[0251]

- 5 x 70% CTX0E03 cultures
- Treat with +/-4OHT, IFN γ , TNF α and TGF β (all at 10ng for 24hrs)
- Detach cells and fix overnight in 2.5% Gluteraldehyde in 0.1M Cacodylate pH7.4
- Cells spun down 300g
- Buffered osmium 2%, 1.5hrs
- Spin, wash water, overnight
- Uranium acetate 2%, 2hrs
- Spin, wash water, 30mins
- Ethanol gradient 20, 35, 50, 70, 80, 90, 100%, over weekend.
- 100% propylene oxide (PO), 1hr
- Spin, 50% Agar LV resin in PO, 1hr
- 75% LV resin/PO 5hrs
- 100% resin overnight at 60°C
- Cool to RT before cutting (60-80nm), Imaged TEM at 200Kv.

Results

[0252] Figure 4A-E shows the electron micrographs of the multivesicular bodies (MVBs) containing exosomes of approximately 30nm - 50nm in diameter. Figure 4F shows microvesicles >100nm in diameter.

Example 8: Production of neural stem cell microparticles from a neural stem cell line.Method

[0253] 5 Sub-confluent flasks containing the same culture of CTX0E03 cells were individually treated with either 10ng/ml TGF- β , 10ng/ml IFN γ , or 10ng/ml TNF α alongside full growth media controls with or without the addition of 4OHT. 72 hours after treatment, the cells were collected using trypzean/EDTA, washed and fixed overnight in 2.5% Gluteraldehyde in 0.1M Cacodylate pH7.4 ready for electron microscopy evaluation.

Results

[0254] The frequency of the occurrence of multivesicular bodies (MVBs) was observed to be altered by the presence of TGF- β , IFN- γ or TNF- α . The frequency was highest in the presence of TGF- β , followed by IFN- γ , followed by TNF- α .

Conclusion

[0255] The production of microparticles from neural stem cells can be stimulated by the addition of the factors TGF- β , IFN- γ or TNF- α . This has the potential for more efficient production of microparticles.

Example 9: Purification, quantification and characterisation of neural stem cell microparticles.Method

[0256] An outline protocol for producing large quantities of microparticles is provided in Figure 5. The main steps are purification, quantification, characterisation, efficacy testing and manufacture.

(1) Purification

Microparticles can be purified from stem cell-conditioned medium by ultracentrifugation, e.g. at 100000 x g for 1-2 hours. Alternative or additional methods for purification of may be used, such as antibody-based methods, e.g. immunoprecipitation, magnetic bead purification, resin-based purification, using specific antibodies.

(2) Quantification

Purified microparticles can be quantified by quantification of total nucleic acid or protein levels, e.g. various PCR or colorimetric protein quantification methods such as such as the BCA assay. Other quantification techniques may alternatively be used, including an electron microscopy grid or an immune-assay using antibodies or antibody fragments that specifically bind to microparticle-specific markers (e.g. ELISA, immunoblotting).

(3) Characterisation

The microparticles can be functionally or structurally characterised. RNA/mRNA/miRNA and protein profiling can be used using methods well known in the art (SDS-PAGE, mass spectrometry, PCR). Constitutively secreted microparticles can be tested and compared to microparticles that have been induced by addition of an inducing agent such as transforming growth factor-beta (TGF- β), interferon-gamma (INF-y) and/or tumour necrosis factor-alpha (TNF- α).

(4) Therapeutic Efficacy

The efficacy of the microparticles can be tested by *in vitro* and *in vivo* assays. For *in vitro* evaluation, neural stem cell microparticles can be added to cultures of monocytes, PBMCs, endothelial cells and/or fibroblasts and the effect of the microparticles on these cells evaluated. Administration of neural stem cell microparticles to suitable animal models can be used to evaluate the *in vivo* efficacy. Clinical trials can be performed to evaluate safety and outcome of neural stem cell microparticles in human subjects.

(5) Manufacture/Scale-Up

Bioreactors, such as the Integra disposable T1000, can be used for the large-scale manufacture of neural stem cell microparticles. The purified microparticles are then formulated as a therapeutic product.

Example 10: Integra CELLINE - Disposable Bioreactor for the production of Micro particles from CTX0E03 cells.

[0257] Efficient micro particle production and harvest from a cell line relies upon maintaining optimal culture conditions for the greatest density of cells. Any restriction in the oxygen or nutrients supplied to the cells or an accumulation of waste metabolic products will limit the life span of the culture, and hence the micro particle production.

[0258] The two-compartment CELLine AD 1000 is designed to accommodate adherent cells attached to a matrix inlay within a small cell compartment, separated from a larger media reservoir by means of a 10kDa semi-permeable mem-

brane. This membrane allows a continuous diffusion of nutrients and removal of waste products, while concentrating any micro particles produced by the cell within the smaller cell compartment. Due to the large volume capacity (1 litre) of the media compartment, the system has the potential to maintain high density cultures for longer periods of time without the need for a media change. The production of exosomes from mesothelioma tumour cell cultures is described in Mitchell *et al*, 2008.

Method

[0259] In order to obtain optimal performance of the CELLline AD1000, place 25ml of complete growth medium (RMM with growth factors and 4OHT) into the medium compartment of the flask to pre-wet the semi-permeable membrane. Allow the flask to sit for 5 minutes at room temperature before coating the matrix inlay with mouse Laminin by adding 15ml of laminin solution (20 μ g/ml in DMEM/F12) to the cell compartment for a minimum of 1 hour at 37°C. Remove the laminin solution and add 15ml of warm DMEM/F12 to the cell compartment to remove any excess laminin. Avoiding the matrix inlay drying, slowly introduce approximately 15x10⁶ CTX0E03 cells in a total of 15ml of complete growth medium. Take care to remove any air bubbles from the cell compartment. Carefully add a further 460ml of complete growth medium to the cell compartment before incubating the flask overnight in 5% CO₂ at 37°C. The next day remove the medium from the cell compartment and replace with 15ml of pre warmed growth medium.

[0260] Every 7 days harvest the microparticles/medium from the cell compartment. Centrifuge the medium at 1500rpm for 5 minutes to remove any cell debris and store at -80°C. Carefully add another 15ml of pre-warmed complete growth medium in to the cell compartment and 485ml of complete growth medium to the medium compartment and incubate for another 7 days. Microparticles were isolated by 100K MWCO filtration. Repeat as necessary.

[0261] Marker characterisations indicated that both purified populations (microvesicles and exosomes) express CD63 and CD81 (determined by FACS - Figure 6). Only the exosomes express the endosomal marker Alix (determined by Western blot, data not shown).

[0262] Figure 8A shows the amount of protein extracted from 15ml of media containing microparticles purified using the Integra system compared to normal culture conditions (3 days T175). Milligrams of protein measured by BCA assay. Figure 8B shows the corresponding quantity of isolated total RNA measured at 260/280nm.

Example 11: Size distribution of Microparticles

[0263] NanoSight analysis was undertaken to determine the particle size and concentration of microvesicles ("mv1" to "mv6") and exosomes ("exo1" to "exo6") isolated from CTX0E03 cells cultured in the Integra Celline system for 1, 2, 3, 4, 5 and 6 weeks. All results are based on 5 replicate measurements.

[0264] Particle size distribution was measured using Nanoparticle Tracking Analysis (NTA). NTA detects the movement of particles in solution and relates it to particle size. Mode and median particle size was calculated for all samples. Exosome samples were analysed using the most sensitive camera settings in order to capture the smallest vesicles. Microvesicle samples were analysed using less sensitive camera settings to prevent over exposure of the larger vesicles. As a result, some smaller vesicles were not detected in the samples. Although smaller vesicles were present in the MV samples, these represent a small percentage of the sample in terms of mass.

[0265] A proportion of Exo1 was labelled with a fluorescent membrane-specific dye (CellMask™) and a combination of NTA analysis with the CellMask™ labelling confirmed that the events detected by NTA correspond to membrane vesicles (data not shown).

[0266] The results are shown in Table 3 below, and in Figure 7.

[0267] The exosomes show a drop in size at week six, from a mode of approximately 110nm to approximately 70nm, or from a median of approximately 130nm to approximately 75nm. The overall size range, from 70nm to 150nm, is consistent with the size of exosomes from other cell types, described in the art. The observed reduction in size of the exosomes to around 70nm diameter after culturing the cells for 6 weeks correlates with the increased efficacy of exosomes isolated from CTX0E03 cells that have been cultured in a multi-compartment bioreactor for 6 weeks correlates, as reported in Example 2 and Figure 3.

[0268] The microvesicles are, as expected, larger, with a mode diameter of approximately 150nm - 200nm, or a median diameter of approximately 180nm - 350nm.

Table 3: Size distribution of CTX0E03 microvesicles and exosomes.

Sample	Count	Dilution	Concentration x10 ¹² /ml	Mode (nm)	Median (nm)
Exo1 (1)	5.204	10000	32.26	107	151
Exo1 (2)	1.734	10000	10.75	135	164

EP 3 033 418 B1

(continued)

Sample	Count	Dilution	Concentration x10¹²/ml	Mode (nm)	Median (nm)
Exo1 (3)	6.55	10000	40.61	108	128
Exo2	14.33	10000	88.85	118	153
Exo3 (1)*	2.52	10000	15.62	89	115
Exo3 (2)	10.06	10000	62.37	115	146
Exo3 (3)	8.98	10000	55.68	128	147
Exo4 (1)	3.04	10000	18.85	111	136
Exo4 (2)	2.89	10000	17.92	110	120
Exo4 (3)	2.77	10000	17.17	116	134
Exo5 (1)	2.34	100	0.15	99	117
Exo5 (2)	2.02	100	0.13	102	124
Exo 5 (3)	2.08	100	0.13	116	127
Exo6 (1)	1.45	100	0.09	68	74
Exo6 (2)	1.19	100	0.07	69	75
MV1 (1)	9.314	200	1.15	183	212
MV1 (2)	10.76	200	1.33	161	214
MV1 (3)	10.738	200	1.33	173	198
MV2	5.89	1000	3.65	177	194
MV3 (1)*	5.68	2000	7.04	150	186
MV3 (2)	11.5	2000	14.26	221	351
MV3 (3)	9.57	2000	11.87	214	270
MV4 (1)	4.894	400	1.21	209	240
MV4 (2)	2.934	1000	1.82	195	212
MV4 (3)	2.55	1000	1.58	184	221
MV5 (1)	1.086	200	0.13	164	237
MV5 (2)	1.458	200	0.18	205	205
MV 5 (3)	1.3	200	0.16	219	210
MV6 (1)	0.346	200	0.04	171	186
MV6 (2)	0.37	200	0.05	168	212

(continued)

Sample	Count	Dilution	Concentration x10 ¹² /ml	Mode (nm)	Median (nm)
Media	0.14	10	0.00	100	149
* large aggregates.					

Example 12: miRNA characterization in CTX0E03 microparticlesMethods**[0269]**

- 3 conditions: CTX0E03 cells in standard culture; microparticles obtained from CTX0E03 cells in standard culture; and purified exosomes derived from CTX0E03 cells in Integra CELLline system (see Examples 10 to 16)
- Investigation of miRNA array using qRT-PCR panel (Qiagen) according to manufacturer's instruction. This assay provides high precision and high sensitivity, with data normalization sensitive to method/choice of reference genes. It does not provide genome wide sequencing.

Results:**[0270]**

A) List of miRNAs with a cp ≤ 35 found in (i) standard CTX0E03 cells, (ii) filtered conditioned medium (0.02-0.2 μ m filter) i.e. microparticles and (iii) exosomes derived from Integra CELLline system (preliminary miRNA qRT-PCR miscript array (Qiagen) results).

B) Arithmetic and geometric mean of the reference (housekeeping) genes

A			
Mature miRNA	CTX0E03 std culture	CM microparticles	CM exosome Integra
hsa-miR-21-5p	19.52	20.9	20.72
hsa-let-7a-5p	22.64	23.11	22.36
hsa-miR-125b-5p	21.64	23.25	21.74
hsa-miR-9-5p	22.58	23.64	22.94
hsa-miR-92a-3p	23.2	23.94	24.01
hsa-miR-24-3p	23.73	24.24	23.83
hsa-miR-20a-5p	23.45	24.43	25.06
hsa-miR-16-5p	23.14	24.72	24.32
hsa-miR-100-5p	23.28	24.74	23.04
hsa-let-7b-5p	24.67	24.75	23.7
hsa-let-7f-5p	23.93	25.09	23.86
hsa-miR-17-5p	24.56	25.24	26.13
hsa-miR-23b-3p	24.3	25.3	24.13
hsa-miR-106b-5p	24.4	25.41	26.16
hsa-miR-222-3p	23.25	25.49	23.17

EP 3 033 418 B1

(continued)

	Mature miRNA	CTX0E03 std culture	CM microparticles	CM exosome Integra
5	hsa-let-7e-5p	24.57	25.58	24.16
	hsa-miR-26a-5p	23.4	25.63	24.2
	hsa-miR-181a-5p	25.16	25.7	24.32
	hsa-miR-125a-5p	23.56	25.75	24.88
10	hsa-miR-103a-3p	24.65	25.8	25.77
	hsa-let-7i-5p	24.37	25.98	24.23
	hsa-miR-99a-5p	24.44	26.05	23.44
15	hsa-let-7c	25.76	26.12	24.07
	hsa-let-7g	25.2	26.15	25.17
	hsa-miR-195-5p	24.72	26.34	25.67
	hsa-miR-93-5p	25.15	26.48	26.06
20	hsa-miR-22-3p	25.03	26.49	25.66
	hsa-miR-20b-5p	26.03	26.86	27.42
	hsa-miR-18a-5p	26.71	26.87	29.06
25	hsa-miR-15b-5p	25.1	26.92	26.43
	hsa-let-7d-5p	26.84	26.96	26.52
	hsa-miR-424-5p	25.56	27.72	26.66
	hsa-miR-15a-5p	26.88	27.89	29.3
30	hsa-miR-130a-3p	27.23	28.26	28.49
	hsa-miR-33a-5p	30.34	28.54	34.18
	hsa-miR-128-	26.94	28.64	27.66
35	hsa-miR-218-5p	27.79	28.68	28.03
	hsa-miR-301a-3p	29.53	28.69	31.57
	hsa-miR-134	28.3	28.76	28.76
	hsa-miR-101-3p	28.44	28.82	31.64
40	hsa-miR-7-5p	29.71	28.82	30.22
	hsa-miR-18b-5p	28.83	28.85	35.47
	hsa-miR-185-5p	28.34	28.99	28.13
45	hsa-miR-378-3p	29.76	29.25	28.97
	hsa-miR-132-3p	28.65	29.32	27.72
	hsa-miR-345-5p	28.49	29.52	29.66
	hsa-miR-219-5p	30.58	29.52	32.7
50	hsa-miR-127-5p	30.05	29.95	31.11
	hsa-miR-146b-5p	30.53	30.54	28.07
	hse-miR-10e-5p	27.1	30.69	28.32
55	hsa-miR-210	29.85	30.83	30.65
	hsa-miR-129-5p	32.51	30.98	31.69
	hsa-miR-137	31.46	31.13	30.95

EP 3 033 418 B1

(continued)

Mature miRNA	CTX0E03 std culture	CM microparticles	CM exosome Integra
hsa-miR-182-5p	28.34	31.64	31.27
hsa-miR-124-3p	33.38	31.71	33.07
hsa-miR-96-5p	29.77	32.27	34.67
hsa-miR-192-5p	31.42	32.42	32.52
hsa-miR-126-3p	31.73	32.44	32.05
hsa-miR-194-5p	31.11	32.49	31.72
hsa-miR-375	33.77	32.94	30.94
hsa-miR-205-5p	35	33.01	32.72
hsa-miR-183-5p	29.88	33.21	31.74
hsa-miR-10b-5p	29.6	33.22	30.79
hsa-miR-302a-3p	29.67	33.6	31.69
hsa-miR-214-3p	34.19	33.76	32.11
hsa-miR-141-3p	35	33.96	34.51
hsa-miR-302c-3p	31.6	34.29	33.93
hsa-miR-196a-5p	35	34.65	35.75
hsa-miR-150-5p	34.59	34.76	34.59
hsa-miR-155-p	32.04	35.75	32.76

B

	CTX0E03 std culture	CM microparticles	CM exosome Integra
Avg. of Arithmetic Mean	23.54	23.82	24.79
Avg. of Geometric Mean	23.48	23.8	24.62

Example 13: CTX0E03 conditioned medium analysis using a protein dot blot

Methods

[0271]

- Conditioned 24hr and 72 hrs conditioned medium (RMM and ITS medium)
- The collected media has been 'concentrated' by dialysis and the proteins biotinylated (typical total protein concentration appears to be 0.5 mg/ml). The media is then incubated with the Raybiotech L507 human protein arrays (total protein concentration 0.1 mg/ml). Following washing and incubation of the array with HRP-conjugated streptavidin, the presence of proteins is detected by chemiluminescence. The array provides qualitative data (i.e. the protein is present, but no indication of its level of expression compared to other proteins).

Results

Cytokine Name	Cytokine Full Name	Function
EDA-A2	ectodysplasin-A2	May be involved in proper formation of skin appendages

(continued)

Cytokine Name	Cytokine Full Name	Function
Galectin-3*	Galectin-3	Galactose-specific lectin which binds IgE. May mediate with the alpha-3, beta-1 integrin the stimulation by CSPG4 of endothelial cells migration.
IGFBP-2	Insulin-like growth factor binding proteins 2	IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture.
IGFBP-rp1/IGFBP-7	Insulin-like Growth Factor Binding Protein Related Protein-1 Insulin-like Growth Factor Binding Protein-7	soluble proteins that bind IGFs with high affinity.
IL-1a†	Interleukin 1 alpha	potent mediator of inflammation and immunity
LECT2†	Leukocyte cell-derived chemotaxin-2	Has a neutrophil chemotactic activity. Also a positive regulator of chondrocyte proliferation.
MCP-1 †	Monocyte chemoattractant protein 1	plays a role in the recruitment of monocytes to sites of injury and infection.
SPARC*	Secreted Protein, Acidic Cysteine-rich-related modular calcium-binding protein 1 [Precursor]	matricellular protein that modulates cell adhesion and proliferation and is thought to function in tissue remodeling and angiogenesis
TIMP-1*	Tissue inhibitor of metalloproteinases-2	Complexes with metalloproteinases (such as collagenases) and irreversibly inactivates them. Also mediates erythropoiesis in vitro; but, unlike IL-3, it is species-specific, stimulating the growth and differentiation of only human and murine erythroid progenitors.
Thrombospondin-1*	Thrombospondin-1	multimodular secreted protein that associates with the extracellular matrix and possesses a variety of biologic functions, including a potent angiogenic activity.
VEGF*	Vascular endothelial growth factor	Growth factor active in angiogenesis, vasculogenesis and endothelial cell growth.
These proteins show expression in some instances -though may also be present in media.		
EGF R/ErbB1	Epidermal growth factor receptor	Receptor for EGF, but also for other members of the EGF family, as TGF-alpha, amphiregulin, betacellulin, heparin-binding EGF-like growth factor
MDC*	A disintegrin and metalloproteinase domain 11 Metalloproteinase-like, disintegrin-like, and cysteine-rich protein MDC	Probable ligand for integrin in the brain. This is a non catalytic metalloprotease-like protein.
Endostatin*	Endostatin	Angiogenesis inhibitor; inhibits endothelial cell migration but may not effect proliferation. May work in balance with VEGF to maintain level of angiogenesis.
Follistatin	Follistatin	Regulates stem cell renewal versus differentiation by inhibiting pro-differentiation proteins

(continued)

Cytokine Name	Cytokine Full Name	Function
Csk†	cytoplasmic tyrosine kinase	Activity is required for interleukin 6 (IL-6) induced differentiation. May play a role in the growth and differentiation of hematopoietic cells. May be involved in signal transduction in endocardial and arterial endothelial cells.
* = angiogenesis † = inflammation		

Example 14: Production of exosomes using the Integra CELLine system.

[0272] CTX0E03 cells were cultured using the Integra CELLine system and exosomes were purified as described in Example 10. The concentration of exosomes purified from the medium using the CELLine system at the 3 week time point, and as a control a standard T175 system as routinely used in the art, was quantified (using a BCA assay to estimate protein content). Figure 9 shows that the production of exosomes using the Integra CELLine system is increased several fold, compared to using conventional culture (T175 flasks).

[0273] Using the Integra CELLine system, CTX0E03 cells were cultured over a 3-week period and medium was harvested at week 1, 2 and 3 for purification and quantification of exosomes, as described in Example 10. Figure 10A shows that the production of microparticles increases exponentially over the 3-week culture period, enabling efficient and large-scale production of microparticles. The concentration of exosomes harvested from a single Integra CELLine flask was then monitored over 1-6 weeks of continuous CTX0E03 culture, with the results shown below and depicted in Figure 10B:

Integra time point	Total quantity of exosomes (ug)	Exosomes ug/ml
Week 1	12	0.80
Week 2	112	7.47
Week 3	88	5.87
Week 4	148	9.87
Week 5	240	16.00
Week 6	440	29.33

[0274] These results show that exosome production is surprisingly enhanced when stem cells are cultured in a multi-compartment bioreactor for weeks, typically at least three weeks.

Example 15: Characterisation of phenotype of cells obtained from the Integra CELLine and the standard (T175) culture system.

[0275] CTX0E03 cells were cultured using the Integra CELLine bioreactor and standard culture, as described in Example 10. Expression of DCX and GFAP protein markers was confirmed using marker-specific antibodies and fluorescence microscopy.

[0276] Expression of DCX, GALC, GFAP, TUBB3, GDNF and IDO markers was detected by qRT-PCR in samples obtained from the cells. Marker expression was compared between microparticles obtained from standard (T175) culture and exosomes obtained from the 3 week cultured Integra CELLine system, assessed against a baseline of the expression level in CTX0E03 cells in standard (T175) culture.

[0277] The inventors observed a striking difference in marker expression of cells obtained from the Integra CELLine system as compared to control cells obtained from standard. Markers of partially-differentiated cells were increased several fold in cells cultured in the Integra CELLine system, compared to control cells obtained from standard cultures (Figure 11). Particularly striking changes are increased expression of the markers DCX1 (doublecortin - a marker for entry into the neural lineage), GFAP (glial fibrillary acidic protein - a marker for entry into the astrocytic lineage), GDNF (glial cell-derived neurotrophic factor) and IDO (indoleamine 2,3-dioxygenase). This indicates that in neural stem cells

cultured in a two-compartment bioreactor partially differentiate into cells of neural (DCX+) or astrocytic (GFAP+) lineage. The expression of DCX and GFAP in the Integra-cultured cells was confirmed by fluorescence microscopy, demonstrating that CTX0E03 cells cultured using the Integra CELLine bioreactor have a more differentiated neuronal phenotype than standard CTX0E03 cells.

Example 16: Characterisation of miRNA expression profiles of exosomes obtained from Integra CELLine cultures and microparticles obtained from standard (T175) cultures.

[0278] CTX0E03 cells were cultured for three weeks using the Integra CELLine culture and in the standard culture in single-compartment T-175 flasks. Exosomes were purified from the Integra culture and microparticles were purified from the standard T-175 culture as described in Example 10. The relative expression levels of various miRNAs expressed in the exosomes and microparticles obtained from either the standard culture or the Integra CELLine system were determined with an miRNA array using qRT-PCR panel (Qiagen) according to manufacturer's instruction, and converted into fold up and down regulation levels as compared to a standard CTX0E03 cell line control group (see Table 4 and Figure 12). These data show a differential miRNA expression profile between exosomes obtained from the Integra CELLine culture system for 3 weeks, microparticles, and cells obtained from the standard single-flask culture.

Table 4: Fold-regulation of miRNAs in microparticles obtained from standard culture or exosomes from the Integra CELLine system, relative to control (CTX0E03 cells).

	Standard Culture (microparticles)	Integra (exosomes)
miRNA	Fold regulation relative to control (CTX0E03 cells)	
hsa-miR-146b-5p	-1.0222	10.5805
hsa-let-7c	-1.6954	4.7678
hsa-miR-99a-5p	-3.5349	3.3714
hsa-miR-132-3p	-1.9163	3.088
hsa-miR-378-3p	1.2731	3.0175
hsa-miR-181a-5p	-1.7431	2.9147
hsa-let-7b-5p	-1.4658	2.7574
hsa-miR-100-5p	-3.208	1.977
hsa-let-7e-5p	-2.7101	1.9274
hsa-miR-23b-3p	-2.3322	1.8834
hsa-miR-185-5p	-1.9119	1.8532
hsa-let-7i-5p	-3.5677	1.8404
hsa-let-7a-5p	-1.851	1.7736
hsa-let-7d-5p	-1.5	1.7654
hsa-let-7g-5p	-2.2527	1.7092
hsa-miR-222-3p	-5.8092	1.6779
hsa-let-7f-5p	-2.8712	1.5948
hsa-miR-218-5p	-1.9611	1.5619
hsa-miR-24-3p	-1.6721	1.5511
hsa-miR-9-5p	-2.2475	1.4109
hsa-miR-126-3p	-2.1263	1.203
hsa-miR-134	-1.6567	1.1783
hsa-miR-128	-3.5842	1.0743
hsa-miR-155-5p	-8.8458	1.0425
hsa-miR-22-3p	-3.4782	-1.0023

EP 3 033 418 B1

(continued)

	Standard Culture (microparticles)	Integra (exosomes)
miRNA	Fold regulation relative to control (CTX0E03 cells)	
hsa-miR-26a-5p	-5.3579	-1.0187
hsa-miR-210	-2.3107	-1.0449
hsa-miR-92a-3p	-1.9885	-1.0693
hsa-miR-93-5p	-3.056	-1.1701
hsa-miR-424-5p	-4.9189	-1.2086
hsa-miR-195-5p	-3.8951	-1.2541
hsa-miR-127-5p	-1.1316	-1.2953
hsa-miR-21-5p	-2.8845	-1.3044
hsa-miR-103a-3p	-2.6482	-1.3287
hsa-miR-16-5p	-3.5267	-1.3692
hsa-miR-125a-5p	-5.1159	-1.434
hsa-miR-10a-5p	-14.4701	-1.434
hsa-miR-10b-5p	-15.1194	-1.4373
hsa-miR-345-5p	-2.5521	-1.4406
hsa-miR-130a-3p	-2.6178	-1.5728
hsa-miR-15b-5p	-4.4025	-1.6058
hsa-miR-20b	-2.1312	-1.6096
hsa-miR-20a-5p	-2.3107	-1.8319
hsa-miR-17-5p	-1.9296	-1.8319
hsa-miR-7-5p	-1.5105	-2.042
hsa-miR-106b-5p	-2.4708	-2.1287
hsa-miR-101-3p	1.4794	-2.4453
hsa-miR-302a-3p	-18.0634	-2.4623
hsa-miR-301a-3p	1.4931	-2.5257
hsa-miR-183-5p	-13.9772	-2.5847
hsa-miR-219-5p	1.6994	-2.7321
hsa-miR-18a-5p	-1.4028	-3.2792
hsa-miR-15a-5p	-2.4766	-3.3714
hsa-miR-182-5p	-12.5099	-4.9588
hsa-miR-33a-5p	2.7927	-9.1472
hsa-miR-96-5p	-7.0047	-18.9396
hsa-miR-18b-5p	-1.3519	-49.18

[0279] Values were calculated from raw data using the following equations:

$$\Delta CT (\text{sample/control}) = \text{Average CT (GOI)} - \text{Average CT (HKG)}$$

$$\text{Fold expression (sample/control)} = 2^{-(\text{Average } \Delta CT)}$$

$$\text{Fold change} = \frac{\text{Fold expression (sample)}}{\text{Fold expression (control)}}$$

If (fold change) > 1 then (fold regulation) = (fold change)

If (fold change) < 1 then (fold regulation) = $-\left(\frac{1}{\text{fold change}}\right)$

[0280] Wherein:

CT = cycle threshold

GOI = gene of interest (investigated miRNA)

HKG = housekeeping genes (reference miRNAs used to normalize the data)

Example 17: Total miRNA analysis

[0281] Cells can shuttle RNA into microparticles determined for release into the extracellular space. This allows the conveyance of genetically encoded messages between cells. We here collectively refer to extracellular RNA as 'shuttle RNA'. We aimed to analyze comprehensively non coding RNA species released by CTX0E03 neural stem cells (NSCs) using Next Generation Sequencing.

[0282] Non coding RNAs are divided in two categories (small and long). Small non coding RNA biotypes include ribosomal RNA (rRNA), small nucleolar (snoRNA), small nuclear RNA (snRNA), microRNA (miRNA), miscellaneous other RNA (misc_RNA, e.g. RMRP, vault RNA, metazoa SRP, and RNY), and long non coding RNA biotypes includes long non-coding RNAs (lncRNAs) and large intergenic non-coding RNAs (lincRNAs).

[0283] Here, we characterized shuttle RNAs, including small and long non coding RNAs, released from NSC derived exosomes and microvesicles (MV) and compared with the RNA contents of the producer NSCs.

A) Total RNA contents in cells, exosomes and microvesicles identified by Agilent RNA bioanalyser

[0284] The RNA in both exosomes and microvesicles mainly consists of small RNA species as shown in Fig. 14. The majority of the nucleotides (nt) was ≤ 200 as shown against the molecular ladder.

B) RNA composition

[0285] Small RNA sequencing libraries were generated to investigate the composition of shuttle and cellular RNA by deep sequencing (Next Generation Sequencing). The results are shown in Figure 15.

C) Deep sequencing of CTX0E03 cell, microvesicle and exosome miRNA expression from standard (T175) cultures.

[0286] Deep sequencing is based on the preparation of a cDNA library following by sequencing and provides information regarding the total sequence read out of different miRNAs in the microvesicles and exosomes. These deep sequence data complement the qRT-PCR array data shown above and provide a comprehensive analysis of the miRNA profile of the cells and microparticles. Unlike the qRT-PCR array analysis, deep sequencing is not restricted to identification of sequences present in the probe array and so the sequences to be identified do not need to be known in advance. Deep sequencing also provides direct read-out and the ability to sequence very short sequences. However, deep sequencing is not suitable for detection of transcripts with low expression.

Method

[0287] The presence of a variety of miRNAs in parental cells and their exosomes (80-100 μ m) and microvesicles (100-1000 μ m), purified by differential centrifugation, was identified by deep sequencing, following construction of 1 tagged miRNA library for each sample.

[0288] Additionally, specific primers for highly shuttled miRNAs (e.g. hsa-miR-1246) were designed and used in real-time reverse transcription PCR (qRT-PCR) to trace exosomes/microvesicles following *in vivo* implantation.

[0289] Deep sequencing was performed by GATC Biotech (Germany) and required the preparation of a tagged miRNA library for each samples followed by sequencing, and miRBase scanning:

- Construction of tagged miRNA libraries (22 to 30 nt)
 - Sequencing libraries were generated by ligation of specific RNA adapter to both 3' and 5' ends for each sample followed by reverse transcription, amplification, and purification of smallRNA libraries (size range of contained smallRNA fraction 22 - 30 nt).
- Sequencing on an Illumina HiSeq 2000 (single read)
 - Sequencing was performed using Illumina HiSeq 2000 (single read). Analysis of one pool could include up to 45,000,000 single read, and each read length is up to 50 bases. Sequencing was quality controlled by using FastQ Files (sequences and quality scores).
- Identification of known miRNAs was performed as followed:
 - RNA adapters were trimmed from resulting sequences and raw data cleaned. Raw data were clustered and for each cluster a number of reads was provided. MiRNAs were identified by miRBase scanning (Ssearch).

Results

[0290] Many microvesicle and exosome miRNAs were enriched relative to the cells, indicating that cells specially sort miRNAs for extracellular release. Furthermore, miRNA contents were similar in both exosomes and microvesicles, indicating a common apparatus of selective miRNA uptake in excreted microvesicles. Without wishing to be bound by theory, this may indicate that miRNA content in secreted microvesicles and exosomes can be used as a fingerprint to identify hNSC subtypes.

[0291] The deep sequencing analysis therefore identified a unique set of miRNAs in both hNSC exosomes and microvesicles not previously reported. MiRNA content in excreted vesicles is similar, but showed a preferential miRNA uptake compared with hNSC. These findings could support biological effects mediated by shuttle miRNA not previously described for hNSC.

[0292] The results are detailed in Tables 5 to 10, below. The data are also depicted in Figure 13, which clearly shows the significantly different miRNA profiles present in the microvesicles and exosomes, compared to the cells. In summary, these data show a massive increase in the amount (read counts) of hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532 in microvesicles and exosomes compared to the cells. Large increases are also seen in hsa-miR-4508, hsa-miR-4516, hsa-miR-3676-5p and hsa-miR-4485. Massive decreases are seen in the amounts (read counts) of certain miRNAs, including hsa-let-7a-5p, hsa-miR-92b-3p, hsa-miR-21-5p, hsa-miR-92a-3p, hsa-miR-10a-5p, hsa-100-5p and hsa-99b-5p.

[0293] The presence of each of hsa-miR-1246, hsa-miR-4488, hsa-miR-4492, hsa-miR-4508, hsa-miR-4516 and hsa-miR-4532 in the exosomes was validated by qRT-PCR (data not shown).

[0294] Plotting the deep sequencing results in the exosomes and microvesicles as relative fold change compared to the cells confirms that hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532 are significantly upregulated in the exosomes and microvesicles compared to the cells. This comparison also shows that miRNA hsa-miR-3195 is the miRNA that is most upregulated, in both exosomes and microvesicles. Although the absolute reads of hsa-miR-3195 are in the range of -40 for exosomes and microvesicles, there is no hsa-miR-3195 present in the cells.

[0295] As noted in Example 16 above, miRNA contents in exosomes, microparticles, and parental cells were also tested and validated using *PCR array* analysis. The following miRNAs were found present by qRT-PCR: hsa-let-7g-5p, hsa-miR-101-3p, hsa-miR-10a-5p, hsa-miR-10b-5p, hsa-miR-125b-5p, hsa-miR-128, hsa-miR-130a-3p, hsa-miR-134, hsa-miR-137, hsa-miR-146b-5p, hsa-miR-15a-5p, hsa-miR-15b-5p, hsa-miR-16-5p, hsa-miR-17-5p, hsa-miR-181a-5p, hsa-miR-182-5p, hsa-miR-185-5p, hsa-miR-18b-5p, hsa-miR-192-5p, hsa-miR-194-5p, hsa-miR-195-5p, hsa-miR-20a-5p, hsa-miR-20b-5p, hsa-miR-210, hsa-miR-21-5p, hsa-miR-218-5p, hsa-miR-219-5p, hsa-miR-222-3p, hsa-miR-22-3p, hsa-miR-23b-3p, hsa-miR-24-3p, hsa-miR-26a-5p, hsa-miR-301a-3p, hsa-miR-302a-3p, hsa-miR-302c-3p, hsa-miR-345-5p, hsa-miR-378a-3p, hsa-miR-7-5p, hsa-miR-92a-3p, hsa-miR-93-5p, hsa-miR-9-5p, hsa-miR-96-5p, and hsa-miR-99a-5p.

EP 3 033 418 B1

Table 5: Cells EH

	Cells: CTX0E03 07EH				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-let-7a-5p	UGAGGUAGUAGGUUGUAUAGUU	1	22	75110
	hsa-miR-10a-5p	UACCCUGUAGAUCCGAAUUUGUG	2	23	52927
10	hsa-miR-100-5p	AACCCGUAGAUCCGAACUUGUG	3	22	52451
	hsa-miR-99b-5p	CACCCGUAGAACCGACCUUGCG	4	22	39457
	hsa-miR-486-5p	UCCUGUACUGAGCUGCCCCGAG	5	22	20310
15	hsa-miR-27b-3p	UUCACAGUGGCUAAGUUCUGC	6	21	16900
	hsa-miR-92a-3p	UAUUGCACUUGUCCCGGCCUGU	7	22	14359
	hsa-miR-191-5p	CAACGGAAUCCCAAAGCAGCUG	8	23	12591
	hsa-miR-21-5p	UAGCUUAUCAGACUGAUGUUGA	9	22	11943
20	hsa-miR-98	UGAGGUAGUAAGUUGUAUUGUU	10	22	11760
	hsa-let-7f-5p	UGAGGUAGUAGAUUGUAUAGUU	11	22	10349
	hsa-miR-26a-5p	UUCAAGUAAUCCAGGAUAGGCU	12	22	9900
25	hsa-miR-92b-3p	UAUUGCACUCGUCCCGGCCUCC	13	22	9794
	hsa-miR-127-3p	UCGGAUCCGUCUGAGCUUGGCU	14	22	7064
	hsa-miR-181a-5p	AACAUUCAACGCUGUCGGUGAGU	15	23	6956
	hsa-miR-182-5p	UUUGGCAAUGGUAGAACUCACACU	16	24	5531
30	hsa-let-7c	UGAGGUAGUAGGUUGUAUGGUU	17	22	5103
	hsa-miR-379-5p	UGGUAGACUAUGGAACGUAGG	18	21	4746
	hsa-miR-146b-5p	UGAGAACUGAAUCCAUAGGCU	19	22	4552
35	hsa-miR-21-3p	CAACACCAGUCGAUGGGCUGU	20	21	4089
	hsa-miR-1246	AAUGGAUUUUUGGAGCAGG	21	19	3973
	hsa-let-7i-5p	UGAGGUAGUAGUUUGUGCUGUU	22	22	3015
	hsa-miR-4532	CCCCGGGGAGCCCGGCG	23	17	2847
40	hsa-miR-183-5p	UAUGGCACUGGUAGAAUUCACU	24	22	2695
	hsa-miR-151a-3p	CUAGACUGAAGCUCCUUGAGG	25	21	2681
	hsa-miR-501-3p	AAUGCACCCGGGCAAGGAUUCU	26	22	2649
45	hsa-let-7e-5p	UGAGGUAGGAGGUUGUAUAGUU	27	22	2449
	hsa-let-7b-5p	UGAGGUAGUAGGUUGUGUGGUU	28	22	2435
	hsa-miR-16-5p	UAGCAGCACGUAAAUAUUGGCG	29	22	2173
	hsa-miR-30a-5p	UGUAAACAUCCUCGACUGGAAG	30	22	2001
50	hsa-miR-30d-5p	UGUAAACAUCCCCGACUGGAAG	31	22	1977
	hsa-miR-409-5p	AGGUUACCCGAGCAACUUUGCAU	32	23	1871
	hsa-miR-22-3p	AAGCUGCCAGUUGAAGAACUGU	33	22	1826
55	hsa-miR-4492	GGGGCUGGGCGCGCGCC	34	17	1754
	hsa-miR-125a-5p	UCCCUGAGACCCUUUAACCUGUGA	35	24	1451

EP 3 033 418 B1

(continued)

5	Cells: CTX0E03 07EH				
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
10	hsa-miR-222-3p	AGCUACAUCUGGCUACUGGGU	36	21	1422
	hsa-miR-151a-5p	UCGAGGAGCUCACAGUCUAGU	37	21	1386
	hsa-miR-181b-5p	AACAUUCAUUGCUGUCGGUGGGU	38	23	1382
	hsa-miR-221-5p	ACCUGGCAUACAAUGUAGAUUU	39	22	1363
	hsa-miR-186-5p	CAAAGAAUUCUCCUUUUGGGCU	40	22	1225
15	hsa-miR-423-5p	UGAGGGGCAGAGAGCGAGACUUU	41	23	1080
	hsa-miR-125b-5p	UCCCUGAGACCCUAACUUGUGA	42	22	1002
	hsa-let-7g-5p	UGAGGUAGUAGUUUGUACAGUU	43	22	959
	hsa-miR-500a-3p	AUGCACCUGGGCAAGGAUUCUG	44	22	923
	hsa-miR-30e-5p	UGUAAACAUCCUUGACUGGAAG	45	22	911
20	hsa-miR-27a-3p	UUCACAGUGGCUAAGUUCCGC	46	21	867
	hsa-miR-409-3p	GAAUGUUGCUCGGUGAACCCCU	47	22	865
	hsa-miR-148b-3p	UCAGUGCAUCACAGAACUUUGU	48	22	856
	hsa-miR-125b-1-3p	ACGGGUUAGGCUCUUGGGAGCU	49	22	851
	hsa-miR-410	AAUAUAACACAGAUGGCCUGU	50	21	848
30	hsa-miR-381	UAUACAAGGGCAAGCUCUCUGU	51	22	842
	hsa-miR-99a-5p	AACCCGUAGAUCCGAUCUUGUG	52	22	773
	hsa-let-7d-5p	AGAGGUAGUAGGUUGCAUAGUU	53	22	765
	hsa-miR-148a-3p	UCAGUGCACUACAGAACUUUGU	54	22	702
	hsa-miR-23a-3p	AUCACAUUGCCAGGGAUUUCC	55	21	654
35	hsa-miR-28-3p	CACUAGAUUGUGAGCUCCUGGA	56	22	593
	hsa-miR-423-3p	AGCUCGGUCUGAGGCCCCUCAGU	57	23	557
	hsa-miR-9-5p	UCUUUGGUUAUCUAGCUGUAUGA	58	23	518
	hsa-miR-23b-3p	AUCACAUUGCCAGGGAUUACC	59	21	508
	hsa-miR-941	CACCCGGCUGUGUGCACAUGUGC	60	23	492
45	hsa-miR-4488	AGGGGGCGGGCUCCGGCG	61	18	485
	hsa-miR-103a-3p	AGCAGCAUUGUACAGGGCUAUGA	62	23	459
	hsa-miR-25-3p	CAUUGCACUUGUCUCGGUCUGA	63	22	436
	hsa-miR-889	UUAUAUCGGACAACCAUUGU	64	21	411
	hsa-miR-378a-3p	ACUGGACUUGGAGUCAGAAGG	65	21	410
50	hsa-miR-30c-5p	UGUAAACAUCCUACACUCUCAGC	66	23	378
	hsa-miR-4485	UAACGGCCGCGGUACCCUAA	67	20	358
	hsa-miR-125b-2-3p	UCACAAGUCAGGCUCUUGGGAC	68	22	352
	hsa-miR-671-3p	UCCGGUUCUCAGGGCUCCACC	69	21	350
	hsa-miR-361-5p	UUAUCAGAAUCUCCAGGGGUAC	70	22	337

EP 3 033 418 B1

(continued)

	Cells: CTX0E03 07EH				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-30e-3p	CUUUCAGUCGGAUGUUUACAGC	71	22	294
	hsa-miR-1271-5p	CUUGGCACCUAGCAAGCACUCA	72	22	288
10	hsa-miR-589-5p	UGAGAACCACGUCUGCUCUGAG	73	22	282
	hsa-miR-374a-5p	UUAUAAUACAACCUGAUAAAGUG	74	22	275
	hsa-miR-769-5p	UGAGACCUCUGGGUUCUGAGCU	75	22	263
	hsa-miR-345-5p	GCUGACUCCUAGUCCAGGGCUC	76	22	249
	hsa-miR-30a-3p	CUUUCAGUCGGAUGUUUGCAGC	77	22	236
	hsa-miR-15b-5p	UAGCAGCACAUCAUGGUUUACA	78	22	229
	hsa-miR-221-3p	AGCUACAUUGUCUGCUGGGUUUC	79	23	225
20	hsa-miR-31-5p	AGGCAAGAUGCUGGCAUAGCU	80	21	213
	hsa-miR-342-3p	UCUCACACAGAAAUCGCACCCGU	81	23	205
	hsa-miR-136-3p	CAUCAUCGUCUCAAAUGAGUCU	82	22	203
25	hsa-miR-493-3p	UGAAGGUCUACUGUGUGCCAGG	83	22	192
	hsa-miR-720	UCUCGCUGGGGCCUCCA	84	17	154
	hsa-miR-7-5p	UGGAAGACUAGUGAUUUUGUUGU	85	23	154
	hsa-miR-130b-3p	CAGUGCAAUGAUGAAAGGGCAU	86	22	150
30	hsa-miR-192-5p	CUGACCUAUGAAUUGACAGCC	87	21	138
	hsa-miR-493-5p	UUGUACAUGGUAGGCUUUCAUU	88	22	115
	hsa-miR-204-5p	UUCCCUUUGUCAUCCU AUGCCU	89	22	113
35	hsa-miR-26b-5p	UUCAAGUAAUUCAGGAUAGGU	90	21	107
	hsa-miR-1307-5p	UCGACCGGACCUCGACCGGCU	91	21	105
	hsa-let-7d-3p	CUAUACGACCUGCUGCCUUUCU	92	22	103
40	hsa-miR-340-5p	UUAUAAAGCAAUGAGACUGAUU	93	22	100
	hsa-miR-134	UGUGACUGGUUGACCAGAGGGG	94	22	99
	hsa-miR-432-5p	UCUUGGAGUAGGUCAUUGGGUGG	95	23	97
	hsa-miR-30b-5p	UGUAAACAUCCUACACUCAGCU	96	22	96
45	hsa-miR-320a	AAAAGCUGGGUUGAGAGGGCGA	97	22	95
	hsa-miR-100-3p	CAAGCUUGUAUCUAUAGGUAUG	98	22	94
	hsa-miR-744-5p	UGC GGGGCUAGGGCUAACAGCA	99	22	89
50	hsa-miR-181a-3p	ACCAUCGACCGUUGAUUGUACC	100	22	86
	hsa-miR-34a-5p	UGGCAGUGUCUUAGCUGGUUGU	101	22	85
	hsa-miR-181a-2-3p	ACCACUGACCGUUGACUGUACC	102	22	81
	hsa-miR-190a	UGAU AUGUUUGAUUAUUAGGU	103	22	79
55	hsa-miR-132-3p	UACAGUCUACAGCCAUGGUCG	104	22	78
	hsa-miR-181c-5p	AACAUUCAACCUGUCGGUGAGU	105	22	76

EP 3 033 418 B1

(continued)

Cells: CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-29a-3p	UAGCACCAUCUGAAAUCGGUUA	106	22	75
hsa-miR-301a-3p	CAGUGCAAUAGUAUUGUCAAGC	107	23	75
hsa-miR-411-5p	UAGUAGACCGUAUAGCGUACG	108	21	75
hsa-miR-128	UCACAGUGAACCGGUCUCUUU	109	21	74
hsa-miR-4516	GGGAGAAGGGUCGGGGC	110	17	74
hsa-miR-425-5p	AAUGACACGAUCACUCCCGUUGA	111	23	72
hsa-miR-130b-5p	ACUCUUUCCCGUUGCACUAC	112	21	71
hsa-miR-130a-3p	CAGUGCAAUGUUAAAAGGGCAU	113	22	67
hsa-miR-30d-3p	CUUUCAGUCAGAUGUUUGCUGC	114	22	65
hsa-miR-654-5p	UGGUGGGCCGCAGAACAUUGUC	115	22	65
hsa-miR-93-5p	CAAAGUGCUGUUCGUGCAGGUAG	116	23	65
hsa-miR-487b	AAUCGUACAGGGUCAUCCACUU	117	22	63
hsa-miR-484	UCAGGCUCAGUCCCCUCCCGAU	118	22	62
hsa-miR-24-3p	UGGCUCAGUUCAGCAGGAACAG	119	22	61
hsa-miR-4677-3p	UCUGUGAGACCAAAGAACUACU	120	22	61
hsa-miR-149-5p	UCUGGCUCCGUGUCUUCACUCCC	121	23	56
hsa-miR-197-3p	UUCACCACCUUCUCCACCCAGC	122	22	56
hsa-miR-96-5p	UUUGGCACUAGCACAUUUUUGCU	123	23	56
hsa-miR-1307-3p	ACUCGGCGUGGCGUCGGUCGUG	124	22	55
hsa-miR-34c-5p	AGGCAGUGUAGUUAGCUGAUUGC	125	23	53
hsa-miR-370	GCCUGCUGGGGUGGAACCUGGU	126	22	52
hsa-miR-148b-5p	AAGUUCUGUUUAACACUCAGGC	127	22	51
hsa-miR-335-5p	UCAAGAGCAAUAACGAAAAUGU	128	23	51
hsa-miR-4461	GAUUGAGACUAGUAGGGCUAGGC	129	23	50
hsa-miR-27a-5p	AGGGCUUAGCUGCUUGUGAGCA	130	22	49
hsa-miR-363-3p	AAUUGCACGGUAUCCAUCUGUA	131	22	47
hsa-miR-431-5p	UGUCUUGCAGGCCGUAUGCA	132	21	47
hsa-miR-877-5p	GUAGAGGAGAUGGCGCAGGG	133	20	46
hsa-miR-550a-5p	AGUGCCUGAGGGAGUAAGAGCCC	134	23	45
hsa-miR-4508	GCGGGGCGUGGGCGCGCG	135	17	44
hsa-miR-541-3p	UGGUGGGCACAGAAUCUGGACU	136	22	42
hsa-miR-135b-5p	UAUGGCUUUUCAUCCUAUGUGA	137	23	40
hsa-miR-140-3p	UACCACAGGGUAGAACCACGG	138	21	39
hsa-miR-362-5p	AAUCCUUGGAACCUAGGUGUGAGU	139	24	37
hsa-miR-455-3p	GCAGUCCAUGGGCAUUAACAC	140	21	37

EP 3 033 418 B1

(continued)

	Cells: CTX0E03 07EH				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-758	UUUGUGACCUGGUCCACUAACC	141	22	37
	hsa-miR-101-3p	UACAGUACUGUGAUAAACUGAA	142	21	36
10	hsa-miR-374b-5p	AUAUAAUACAACCUGCUAAGUG	143	22	36
	hsa-miR-148a-5p	AAAGUUCUGAGACACUCCGACU	144	22	35
	hsa-miR-17-5p	CAAAGUGCUUACAGUGCAGGUAG	145	23	35
15	hsa-miR-20a-5p	UAAAGUGCUUAUAGUGCAGGUAG	146	23	35
	hsa-miR-874	CUGCCCUGGCCCGAGGGACCGA	147	22	35
	hsa-miR-193b-3p	AACUGGCCCUCAAAGUCCCGCU	148	22	34
	hsa-miR-548ah-3p	CAAAAACUGCAGUUACUUUUGC	149	22	34
20	hsa-miR-539-3p	AUCAUACAAGGACAAUUUCUUU	150	22	33
	hsa-miR-421	AUCAACAGACAUUAAUUGGGCGC	151	23	31
	hsa-miR-28-5p	AAGGAGCUCACAGUCUAUUGAG	152	22	30
25	hsa-miR-485-3p	GUCAUACACGGCUCUCCUCUCU	153	22	29
	hsa-miR-2467-5p	UGAGGCUCUGUUAGCCUUGGCUC	154	23	26
	hsa-miR-4449	CGUCCCGGGGCUGCGCGAGGCA	155	22	26
	hsa-miR-24-2-5p	UGCCUACUGAGCUGAAACACAG	156	22	25
30	hsa-miR-181d	AACAUUCAUUGUUGUCGGUGGGU	157	23	24
	hsa-miR-323a-3p	CACAUUACACGGUCGACCUCU	158	21	24
	hsa-miR-106b-3p	CCGCACUGUGGGUACUUGCUGC	159	22	23
35	hsa-miR-125a-3p	ACAGGUGAGGUUCUUGGGAGCC	160	22	23
	hsa-miR-330-5p	UCUCUGGGCCUGUGUCUUAGGC	161	22	23
	hsa-miR-1275	GUGGGGGAGAGGCUGUC	162	17	22
	hsa-miR-19b-3p	UGUGCAAAUCCAUGCAAAACUGA	163	23	22
40	hsa-miR-301b	CAGUGCAAUGAUUUGUCAAAAGC	164	23	21
	hsa-miR-485-5p	AGAGGCUGGCCGUGAUGAAUUC	165	22	21
	hsa-miR-29b-3p	UAGCACCAUUUGAAAUCAGUGUU	166	23	20
45	hsa-miR-3158-3p	AAGGGCUUCCUCUCUGCAGGAC	167	22	20
	hsa-miR-431-3p	CAGGUCGUCUUGCAGGGCUUCU	168	22	20
	hsa-miR-454-3p	UAGUGCAAUAUUGCUUAUAGGGU	169	23	20
50	hsa-miR-106b-5p	UAAAGUGCUGACAGUGCAGAU	170	21	19
	hsa-miR-1973	ACCGUGCAAAGGUAGCAUA	171	19	19
	hsa-miR-31-3p	UGCUAUGCCAACAUUUGCCAU	172	22	19
	hsa-miR-374a-3p	CUUAUCAGAUUGUAUUGUAAUU	173	22	19
55	hsa-miR-433	AUCAUGAUGGGCUCCUCGGUGU	174	22	19
	hsa-miR-4417	GGUGGGCUUCCCGGAGGG	175	18	19

EP 3 033 418 B1

(continued)

Cells: CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-143-3p	UGAGAUGAAGCACUGUAGCUC	176	21	18
hsa-miR-19a-3p	UGUGCAAAUCUAUGCAAAACUGA	177	23	18
hsa-miR-532-5p	CAUGCCUUGAGUGUAGGACCGU	178	22	18
hsa-miR-561-5p	AUCAAGGAUCUUAACUUGCC	179	22	18
hsa-miR-663b	GGUGGCCCGGCCGUGCCUGAGG	180	22	18
hsa-miR-1301	UUGCAGCUGCCUGGGAGUGACUUC	181	24	17
hsa-miR-299-3p	UAUGUGGGAUGGUAACCGCUU	182	22	17
hsa-miR-9-3p	AUAAAGCUAGAUAAACCGAAAGU	183	22	17
hsa-miR-17-3p	ACUGCAGUGAAGGCACUUGUAG	184	22	15
hsa-miR-376c	AACAUAGAGGAAAUUCCACGU	185	21	15
hsa-miR-424-5p	CAGCAGCAAUUCAUGUUUUGAA	186	22	15
hsa-miR-660-5p	UACCCAUUGCAUAUCGGAGUUG	187	22	15
hsa-miR-153	UUGCAUAGUCACAAAAGUGAUC	188	22	14
hsa-miR-3605-5p	UGAGGAUGGAUAGCAAGGAAGCC	189	23	14
hsa-miR-3687	CCCGGACAGGCGUUCGUGCGACGU	190	24	14
hsa-miR-4284	GGGCUCACAUCACCCCAU	191	18	14
hsa-miR-455-5p	UAUGUGCCUUUGGACUACAUCG	192	22	14
hsa-miR-543	AAACAUUCGCGGUGCACUUCUU	193	22	14
hsa-miR-1276	UAAAGAGCCCUGUGGAGACA	194	20	13
hsa-miR-330-3p	GCAAAGCACACGGCCUGCAGAGA	195	23	13
hsa-miR-369-3p	AAUAAUACAUGGUUGAUCUUU	196	21	13
hsa-miR-4786-5p	UGAGACCAGGACUGGAUGCACC	197	22	13
hsa-miR-548k	AAAAGUACUUGCGGAUUUUGCU	198	22	13
hsa-miR-1226-3p	UCACCAGCCCUGUGUCCCCUAG	199	22	12
hsa-miR-188-3p	CUCCCAUAUGCAGGGUUUGCA	200	21	12
hsa-miR-27b-5p	AGAGCUUAGCUGAUUGGUGAAC	201	22	12
hsa-miR-377-5p	AGAGGUUGCCCUUGGUGAAUUC	202	22	12
hsa-miR-487a	AAUCAUACAGGGACAUCAGUU	203	22	12
hsa-miR-92a-1-5p	AGGUUGGGAUCGGUUGCAAUGCU	204	23	12
hsa-miR-135b-3p	AUGUAGGGCUAAAAGCCAUGGG	205	22	11
hsa-miR-218-5p	UUGUGCUUGAUCUAACCAUGU	206	21	11
hsa-miR-3943	UAGCCCCCAGGCUUCACUUGGCG	207	23	11
hsa-miR-92b-5p	AGGGACGGGACGCGGUGCAGUG	208	22	11
hsa-miR-1185-1-3p	AUAUACAGGGGGAGACUCUUAU	209	22	10
hsa-miR-1273g-3p	ACCACUGCACUCCAGCCUGAG	210	21	10

EP 3 033 418 B1

(continued)

5	Cells: CTX0E03 07EH			
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH
10	hsa-miR-2355-5p	AUCCCCAGAUACAAUGGACAA	211	21
	hsa-miR-23a-5p	GGGGUUCUGGGGAUGGGAUUU	212	22
15	hsa-miR-30c-1-3p	CUGGGAGAGGGUUGUUUACUCC	213	22
	hsa-miR-329	AACACACCUGGUUAACCUCUUU	214	22
20	hsa-miR-337-3p	CUCCUAUAUGAUGCCUUUCUUC	215	22
	hsa-miR-3609	CAAAGUGAUGAGUAAUACUGGCUG	216	24
25	hsa-miR-378a-5p	CUCCUGACUCCAGGUCCUGUGU	217	22
	hsa-miR-3929	GAGGCUGAUGUGAGUAGACCACU	218	23
30	hsa-miR-4745-5p	UGAGUGGGGCUCCCGGGACGGCG	219	23
	hsa-miR-5096	GUUUCACCAUGUUGGUCAGGC	220	21
35	hsa-miR-656	AAUAUUUAACAGUCAACCUCU	221	21
	hsa-let-7a-3p	CUAUACAAUCUACUGUCUUUC	222	21
40	hsa-miR-15a-5p	UAGCAGCACAAUUGGUUUGUG	223	22
	hsa-miR-185-5p	UGGAGAGAAAGGCAGUCCUGA	224	22
45	hsa-miR-25-5p	AGGCGGAGACUUGGGCAAUUG	225	21
	hsa-miR-3065-5p	UCAACAAAAUCACUGAUGCUGGA	226	23
50	hsa-miR-3176	ACUGGCCUGGGACUACCGG	227	19
	hsa-miR-339-3p	UGAGCGCCUCGACGACAGAGCCG	228	23
55	hsa-miR-374b-3p	CUUAGCAGGUUGUAUUAUCAUU	229	22
	hsa-miR-4435	AUGGCCAGAGCUCACACAGAGG	230	22
60	hsa-miR-4448	GGCUCCUUGGUCUAGGGGUA	231	20
	hsa-miR-4497	CUCCGGGACGGCUGGGC	232	17
65	hsa-miR-4521	GCUAAGGAAGUCCUGUGCUCAG	233	22
	hsa-miR-539-5p	GGAGAAUUUAUCCUUGGUGUGU	234	22
70	hsa-miR-548ah-5p	AAAAGUGAUUGCAGUGUUUG	235	20
	hsa-miR-1910	CCAGUCCUGUGCCUGCCGCCU	236	21
75	hsa-miR-376a-3p	AUCAUAGAGGAAAAUCCACGU	237	21
	hsa-miR-382-5p	GAAGUUGUUCGUGGUGGAUUCG	238	22
80	hsa-miR-3940-3p	CAGCCCGGAUCCCAGCCCACUU	239	22
	hsa-miR-494	UGAAACAUACACGGGAAACCUC	240	22
85	hsa-miR-495	AAACAAACAUGGUGCACUUCUU	241	22
	hsa-miR-545-3p	UCAGCAAACAUUUUUGUGUGC	242	22
90	hsa-miR-99b-3p	CAAGCUCGUGUCUGUGGGUCCG	243	22
	hsa-miR-1197	UAGGACACAUGGUCUACUUCU	244	21
95	hsa-miR-181b-3p	CUCACUGAACAAUGAAUGCAA	245	21

EP 3 033 418 B1

(continued)

	Cells: CTX0E03 07EH				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-212-5p	ACCUUGGCUCUAGACUGCUUACU	246	23	7
	hsa-miR-3200-3p	CACCUUGCGCUACUCAGGUCUG	247	22	7
10	hsa-miR-340-3p	UCCGUCUCAGUUACUUUAUAGC	248	22	7
	hsa-miR-3607-5p	GCAUGUGAUGAAGCAAAUCAGU	249	22	7
	hsa-miR-361-3p	UCCCCCAGGUGUGAUUCUGAUUU	250	23	7
15	hsa-miR-3656	GGCGGGUGCGGGGGUGG	251	17	7
	hsa-miR-532-3p	CCUCCCACACCCAAGGCUUGCA	252	22	7
	hsa-miR-574-3p	CACGCUCAUGCACACACCCACA	253	22	7
	hsa-miR-107	AGCAGCAUUGUACAGGGCUAUCA	254	23	6
20	hsa-miR-127-5p	CUGAAGCUCAGAGGGCUCUGAU	255	22	6
	hsa-miR-18a-5p	UAAGGUGCAUCUAGUGCAGAUAG	256	23	6
	hsa-miR-26a-2-3p	CCUAUUCUUGAUUACUUGUUUC	257	22	6
25	hsa-miR-296-5p	AGGGCCCCCCCCUCAAUCCUGU	258	21	6
	hsa-miR-3648	AGCCGCGGGGAUCGCCGAGGG	259	21	6
	hsa-miR-382-3p	AAUCAUUCACGGACAACACUU	260	21	6
	hsa-miR-3939	UACGCGCAGACCACAGGAUGUC	261	22	6
30	hsa-miR-432-3p	CUGGAUGGCUCCUCCAUGUCU	262	21	6
	hsa-miR-4423-5p	AGUUGCCUUUUUGUUCCCAUGC	263	22	6
	hsa-miR-4466	GGGUGCGGGCCGGCGGGG	264	18	6
35	hsa-miR-454-5p	ACCCUAUCAAUUUGUCUCUGC	265	22	6
	hsa-miR-4746-5p	CCGGUCCCAGGAGAACCUGCAGA	266	23	6
	hsa-miR-496	UGAGUAUUACAUGGCCAAUCUC	267	22	6
	hsa-miR-548o-3p	CCAAAACUGCAGUUACUUUUGC	268	22	6
40	hsa-miR-1248	ACCUUCUUGUAUAAGCACUGUGCUGAAA	269	27	5
	hsa-miR-1254	AGCCUGGAAGCUGGAGCCUGCAGU	270	24	5
	hsa-miR-1296	UUAGGGCCCUGGCUCCAUCUCC	271	22	5
45	hsa-miR-136-5p	ACUCCAUUUGUUUUGAUGAUGGA	272	23	5
	hsa-miR-199a-5p	CCCAGUGUUCAGACUACCUGUUC	273	23	5
	hsa-miR-296-3p	GAGGGUUGGGUGGAGGCUCUCC	274	22	5
50	hsa-miR-3177-3p	UGCACGGCACUGGGGACACGU	275	21	5
	hsa-miR-324-3p	ACUGCCCCAGGUGCUGCUGG	276	20	5
	hsa-miR-337-5p	GAACGGCUUCAUACAGGAGUU	277	21	5
	hsa-miR-342-5p	AGGGGUGCUAUCUGUGAUUGA	278	21	5
55	hsa-miR-365b-3p	UAAUGCCCCUAAAAUCCUUAU	279	22	5
	hsa-miR-3676-5p	AGGAGAUCCUGGGUU	280	15	5

EP 3 033 418 B1

(continued)

5	Cells: CTX0E03 07EH			
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH
10	hsa-miR-502-3p	AAUGCACCUGGGCAAGGAUUA	281	22
	hsa-miR-505-3p	CGUCAACACUUGCUGGUUUCCU	282	22
15	hsa-miR-550a-3p	UGUCUUACUCCCUCAGGCACAU	283	22
	hsa-miR-5587-3p	GCCCCGGGCAGUGUGAUCAUC	284	21
20	hsa-miR-641	AAAGACAUAGGAUAGAGUACCUC	285	24
	hsa-miR-655	AUAAUACAUGGUUAACCUCUUU	286	22
25	hsa-miR-664-3p	UAUUCAUUUAUCCCCAGCCUACA	287	23
	hsa-miR-671-5p	AGGAAGCCCUGGAGGGGCUGGAG	288	23
30	hsa-miR-760	CGGCUCUGGGUCUGUGGGGA	289	20
	hsa-let-7e-3p	CUAUACGGCCUCCUAGCUUUC	290	22
35	hsa-miR-1268a	CGGGCGUGGUGGUGGGGG	291	18
	hsa-miR-1273f	GGAGAUUGGAGGUUGCAGUG	292	19
40	hsa-miR-1286	UGCAGGACCAAGAUGAGCCCU	293	21
	hsa-miR-1291	UGGCCCUGACUGAAGACCAGCAGU	294	24
45	hsa-miR-141-3p	UAACACUGUCUGGUAAAGAUGG	295	22
	hsa-miR-1468	CUCCGUUUGCCUGUUUCGCUG	296	21
50	hsa-miR-328	CUGGCCUCUCUGCCCUUCCGU	297	22
	hsa-miR-424-3p	CAAAACGUGAGGCGCUGCUAU	298	21
55	hsa-miR-4454	GGAUCCGAGUCACGGCACCA	299	20
	hsa-miR-4463	GAGACUGGGUGGGGCC	300	17
60	hsa-miR-4671-3p	UUAGUGCAUAGUCUUUGGUCU	301	21
	hsa-miR-4775	UUAAUUUUUUGUUUCGGUCACU	302	22
65	hsa-miR-500a-5p	UAAUCCUUGCUACCUGGGUGAGA	303	23
	hsa-miR-548b-5p	AAAAGUAAUUGUGUUUUGGCC	304	22
70	hsa-miR-573	CUGAAGUGAUGUGUAACUGAUCAG	305	24
	hsa-miR-576-5p	AUUCUAAUUUCUCCACGUCUUU	306	22
75	hsa-miR-625-3p	GACUAUAGAACUUUCCCCCUCA	307	22
	hsa-miR-652-3p	AAUGGCGCCACUAGGGUUGUG	308	21
80	hsa-miR-665	ACCAGGAGGCUGAGGCCCCU	309	20
	hsa-miR-766-3p	ACUCCAGCCCCACAGCCUCAGC	310	22
85	hsa-miR-935	CCAGUUACCGCUUCCGCUACCGC	311	23
	hsa-miR-937	AUCCGCGCUCUGACUCUCUGCC	312	22
90	hsa-miR-1180	UUUCCGGCUCGCGUGGGUGUGU	313	22
	hsa-miR-1185-2-3p	AUAUACAGGGGAGACUCUCAU	314	22
95	hsa-miR-132-5p	ACCGUGGCUUUCGAUUGUUACU	315	22

EP 3 033 418 B1

(continued)

Cells: CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-16-2-3p	CCAAUAUUACUGUGCUGCUUUA	316	22	3
hsa-miR-20b-5p	CAAAGUGCUCAUAGUGCAGGUAG	317	23	3
hsa-miR-2116-3p	CCUCCCAUGCCAAGAACUCCC	318	21	3
hsa-miR-299-5p	UGGUUUACCGUCCCACAUACAU	319	22	3
hsa-miR-30b-3p	CUGGGAGGUGGAUGUUUACUUC	320	22	3
hsa-miR-30c-2-3p	CUGGGAGAAGGCUGUUUACUCU	321	22	3
hsa-miR-3187-3p	UUGGCCAUGGGGCUGCGCGG	322	20	3
hsa-miR-3615	UCUCUCGGCUCCUCGCGGCUC	323	21	3
hsa-miR-3620	UCACCCUGCAUCCCGCACCCAG	324	22	3
hsa-miR-3654	GACUGGACAAGCUGAGGAA	325	19	3
hsa-miR-3662	GAAAAUGAUGAGUAGUGACUGAUG	326	24	3
hsa-miR-3681-5p	UAGUGGAUGAUGCACUCUGUGC	327	22	3
hsa-miR-4286	ACCCACUCCUGGUACC	328	17	3
hsa-miR-4640-3p	CACCCCCUGUUUCCUGGCCAC	329	22	3
hsa-miR-4717-3p	ACACAUGGGUGGCUGUGGCCU	330	21	3
hsa-miR-542-3p	UGUGACAGAUUGAUAAACUGAAA	331	22	3
hsa-miR-5584-5p	CAGGGAAAUGGGAAGAACUAGA	332	22	3
hsa-miR-570-3p	CGAAAACAGCAAUUACCUUUGC	333	22	3
hsa-miR-574-5p	UGAGUGUGUGUGUGUGAGUGUGU	334	23	3
hsa-miR-628-3p	UCUAGUAAGAGUGGCAGUCGA	335	21	3
hsa-miR-654-3p	UAUGUCUGCUGACCAUCACCUU	336	22	3
hsa-miR-769-3p	CUGGGAUCUCCGGGGUCUUGGUU	337	23	3
hsa-miR-943	CUGACUGUUGCCGUCCUCCAG	338	21	3
hsa-let-7b-3p	CUAUACAACCUACUGCCUCCCC	339	22	2
hsa-miR-1244	AAGUAGUUGGUUUGUAUGAGAUGGUU	340	26	2
hsa-miR-1255a	AGGAUGAGCAAAGAAAGUAGAUU	341	23	2
hsa-miR-1273e	UUGCUUGAACCCAGGAAGUGGA	342	22	2
hsa-miR-1289	UGGAGUCCAGGAUCUGCAUUUU	343	23	2
hsa-miR-152	UCAGUGCAUGACAGAACUUGG	344	21	2
hsa-miR-194-5p	UGUAACAGCAACUCCAUGUGGA	345	22	2
hsa-miR-195-5p	UAGCAGCACAGAAAUUUGGC	346	21	2
hsa-miR-200c-3p	UAAUACUGCCGGGUAAUGAUGGA	347	23	2
hsa-miR-212-3p	UACAGUCUCCAGUCACGGCC	348	21	2
hsa-miR-222-5p	CUCAGUAGCCAGUGUAGAUCU	349	22	2
hsa-miR-3065-3p	UCAGCACCAGGAUAUUGUUGGAG	350	23	2

EP 3 033 418 B1

(continued)

Cells: CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-3115	AUAUGGGUUUACUAGUUGGU	351	20	2
hsa-miR-3126-5p	UGAGGGACAGAUGCCAGAAGCA	352	22	2
hsa-miR-3174	UAGUGAGUUAGAGAUAGCAGAGCC	353	23	2
hsa-miR-324-5p	CGCAUCCCCUAGGGCAUUGGUGU	354	23	2
hsa-miR-33a-5p	GUGCAUUGUAGUUGCAUUGCA	355	21	2
hsa-miR-3677-3p	CUCGUGGGCUCUGGCCACGGCC	356	22	2
hsa-miR-369-5p	AGAUCGACCGUGUUUAUUAUUCGC	357	22	2
hsa-miR-425-3p	AUCGGGAAUGUCGUGUCCGCCC	358	22	2
hsa-miR-4426	GAAGAUGGACGUACUUU	359	17	2
hsa-miR-4467	UGGCGGCGGUAGUUAUGGGCUU	360	22	2
hsa-miR-4482-3p	UUUCUAUUUCUCAGUGGGGCUC	361	22	2
hsa-miR-4515	AGGACUGGACUCCCGGCAGCCC	362	22	2
hsa-miR-4792	CGGUGAGCGCUCGCUGGC	363	18	2
hsa-miR-659-5p	AGGACCUUCCCUGAACCAAGGA	364	22	2
hsa-miR-663a	AGGCGGGGCGCCGCGGGACCGC	365	22	2
hsa-miR-940	AAGGCAGGGCCCCCGCUCCCC	366	21	2
hsa-miR-99a-3p	CAAGCUCGCUUCUAUGGGUCUG	367	22	2
hsa-miR-1185-5p	AGAGGAUACCCUUUGUAUGUU	368	21	1
hsa-miR-1225-3p	UGAGCCCCUGUGCCGCCCCAG	369	22	1
hsa-miR-1237	UCCUUCUGCUCCGUCCCCAG	370	21	1
hsa-miR-1252	AGAAGGAAAU UGAAU UCAU U UA	371	22	1
hsa-miR-1257	AGUGAAUGAUGGGUUCUGACC	372	21	1
hsa-miR-1260b	AUCCCACCACUGCCACCAU	373	19	1
hsa-miR-1273d	GAACCAUGAGGUUGAGGCUGCAGU	374	25	1
hsa-miR-1290	UGGAUUUUUGGAUCAGGGA	375	19	1
hsa-miR-1306-3p	ACGUUGGCUCUGGUGGUG	376	18	1
hsa-miR-1321	CAGGGAGGUGAAUGUGAU	377	18	1
hsa-miR-1343	CUCCUGGGGCCCCGCACUCUCGC	378	22	1
hsa-miR-138-5p	AGCUGGUGUUGUGAAUCAGGCCG	379	23	1
hsa-miR-140-5p	CAGUGGUUUUACCCUAUGGUAG	380	22	1
hsa-miR-146b-3p	UGCCUGUGGACUCAGUUCUGG	381	22	1
hsa-miR-186-3p	GCCCAAAGGUGAAUUUUUUGGG	382	22	1
hsa-miR-1908	CGGCGGGGACGGCGAUUGGUC	383	21	1
hsa-miR-1915-3p	CCCCAGGGCGACGCGGCGGG	384	20	1
hsa-miR-1915-5p	ACCUUGCCUUGCUGCCCGGGCC	385	22	1

EP 3 033 418 B1

(continued)

Cells: CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-193a-3p	AACUGGCCUACAAAGUCCCAGU	386	22	1
hsa-miR-19b-1-5p	AGUUUUGCAGGUUUGCAUCCAGC	387	23	1
hsa-miR-208b	AUAAGACGAACAAAAGGUUUGU	388	22	1
hsa-miR-2110	UUGGGGAAACGGCCGCUGAGUG	389	22	1
hsa-miR-219-1-3p	AGAGUUGAGUCUGGACGUCCCG	390	22	1
hsa-miR-26b-3p	CCUGUUCUCCAUIJACUUGGCUC	391	22	1
hsa-miR-2964a-3p	AGAAUUGCGUUUGGACAAUCAGU	392	23	1
hsa-miR-29a-5p	ACUGAUUUUCUUUUGGUGUUCAG	393	22	1
hsa-miR-3126-3p	CAUCUGGCAUCCGUCACACAGA	394	22	1
hsa-miR-3130-3p	GCUGCACCGGAGACUGGGUAA	395	21	1
hsa-miR-3130-5p	UACCCAGUCUCCGGUGCAGCC	396	21	1
hsa-miR-3140-5p	ACCUGAAUUACCAAAAGCUUU	397	21	1
hsa-miR-3155a	CCAGGCUCUGCAGUGGGAACU	398	21	1
hsa-miR-3157-3p	CUGCCCUAGUCUAGCUGAAGCU	399	22	1
hsa-miR-3180-3p	UGGGGCGGAGCUUCCGGAGGCC	400	22	1
hsa-miR-323b-5p	AGGUUGUCCGUGGUGAGUUCGCA	401	23	1
hsa-miR-339-5p	UCCCUGUCCUCCAGGAGCUCACG	402	23	1
hsa-miR-34a-3p	CAAUCAGCAAGUAUACUGCCCU	403	22	1
hsa-miR-34b-3p	CAAUCACUAACUCCACUGCCAU	404	22	1
hsa-miR-34c-3p	AAUCACUAACCACACGGCCAGG	405	22	1
hsa-miR-3658	UUUAAGAAAACACCAUGGAGAU	406	22	1
hsa-miR-365a-5p	AGGGACUUUUGGGGGCAGAUGUG	407	23	1
hsa-miR-3676-3p	CCGUGUUUCCCCACGCUUU	408	20	1
hsa-miR-3691-5p	AGUGGAUGAUGGAGACUCGGUAC	409	23	1
hsa-miR-376a-5p	GUAGAUUCUCCUUCUAUGAGUA	410	22	1
hsa-miR-378g	ACUGGGCUUGGAGUCAGAAG	411	20	1
hsa-miR-3909	UGUCCUCUAGGGCCUGCAGUCU	412	22	1
hsa-miR-3928	GGAGGAACCUUGGAGCUUCGGC	413	22	1
hsa-miR-3942-3p	UUUCAGAUAACAGUAUUACAU	414	21	1
hsa-miR-3944-5p	UGUGCAGCAGGCCAACCGAGA	415	21	1
hsa-miR-3960	GGCGGCGGCGGAGGCGGGGG	416	20	1
hsa-miR-4326	UGUUCUCUGUCUCCAGAC	417	20	1
hsa-miR-4444	CUCGAGUUGGAAGAGGCG	418	18	1
hsa-miR-4450	UGGGGAUUUGGAGAAGUGGUGA	419	22	1
hsa-miR-4642	AUGGCAUCGUCCCCUGGUGGCU	420	22	1

EP 3 033 418 B1

(continued)

Cells: CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-4668-5p	AGGGAAAAAAAAAAGGAUUUGUC	421	23	1
hsa-miR-4673	UCCAGGCAGGAGCCGGACUGGA	422	22	1
hsa-miR-4688	UAGGGGCAGCAGAGGACCUGGG	423	22	1
hsa-miR-4700-3p	CACAGGACUGACUCCUCACCCCAGUG	424	26	1
hsa-miR-4731-3p	CACACAAGUGGCCCCCAACACU	425	22	1
hsa-miR-4749-3p	CGCCCCUCCUGCCCCACAG	426	20	1
hsa-miR-4769-5p	GGUGGGAUGGAGAGAAGGUAUGAG	427	24	1
hsa-miR-4800-5p	AGUGGACCGAGGAAGGAAGGA	428	21	1
hsa-miR-491-5p	AGUGGGGAACCCUCCAUGAGG	429	22	1
hsa-miR-501-5p	AAUCCUUUGUCCUGGGUGAGA	430	22	1
hsa-miR-5092	AAUCCACGCUGAGCUUGGCAUC	431	22	1
hsa-miR-541-5p	AAAGGAUUCUGCUGUCGGUCCACU	432	25	1
hsa-miR-542-5p	UCGGGGAUCAUCAUGUCACGAGA	433	23	1
hsa-miR-551b-3p	GCGACCAUACUUGGUUUCAG	434	21	1
hsa-miR-5690	UCAGCUACUACCUCUAUUAGG	435	21	1
hsa-miR-577	UAGAUAAAAUUAUUGGUACCUG	436	21	1
hsa-miR-584-3p	UCAGUUCAGGCCAACCCAGGCU	437	22	1
hsa-miR-589-3p	UCAGAACAAAUGCCGGUCCAGAGA	438	24	1
hsa-miR-616-5p	ACUCAAAACCCUUCAGUGACUU	439	22	1
hsa-miR-628-5p	AUGCUGACAUUUUACUAGAGG	440	22	1
hsa-miR-629-5p	UGGGUUUACGUUGGGAGAACU	441	21	1
hsa-miR-644b-3p	UUCAUUUGCCUCCAGCCUACA	442	22	1
hsa-miR-664-5p	ACUGGCUAGGGAAAAUGAUUGGAU	443	24	1
hsa-miR-922	GCAGCAGAGAAUAGGACUACGUC	444	23	1

Table 6: Cells EI

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-let-7a-5p	UGAGGUAGUAGGUUGUAUAGUU	1	22	305060
hsa-miR-92b-3p	UAUUGCACUCGUCCCGGCCUCC	13	22	242715
hsa-miR-21-5p	UAGCUUAUCAGACUGAUUGUA	9	22	154626
hsa-miR-92a-3p	UAUUGCACUUGUCCCGGCCUGU	7	22	137412
hsa-miR-127-3p	UCGGAUCCGUCUGAGCUUGGCU	14	22	110806
hsa-miR-100-5p	AACCCGUAGAUCGGAACUUGUG	3	22	109290

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-27b-3p	UUCACAGUGGCUAAGUUCUGC	6	21	91902
hsa-miR-191-5p	CAACGGAAUCCCAAAAGCAGCUG	8	23	89150
hsa-miR-26a-5p	UUCAAGUAAUCCAGGAUAGGCU	12	22	88724
hsa-miR-99b-5p	CACCCGUAGAACCGACCUUGCG	4	22	87399
hsa-let-7f-5p	UGAGGUAGUAGAUUGUAUAGUU	11	22	78395
hsa-miR-181a-5p	AACAUUCAACGCUGUCGGUGAGU	15	23	47686
hsa-miR-486-5p	UCCUGUACUGAGCUGCCCCGAG	5	22	41639
hsa-miR-30a-5p	UGUAAACAUCCUCGACUGGAAG	30	22	35465
hsa-miR-98	UGAGGUAGUAAGUUGUAUUGUU	10	22	30440
hsa-miR-151a-3p	CUAGACUGAAGCUCCUUGAGG	25	21	29047
hsa-miR-21-3p	CAACACCAGUCGAUGGGCUGU	20	21	27733
hsa-miR-30d-5p	UGUAAACAUCCCCGACUGGAAG	31	22	27307
hsa-let-7c	UGAGGUAGUAGGUUGUAUGGUU	17	22	27224
hsa-miR-10a-5p	UACCCUGUAGAUCCGAAUUUGUG	2	23	26908
hsa-miR-22-3p	AAGCUGCCAGUUGAAGAACUGU	33	22	26456
hsa-miR-182-5p	UUUGGCAAUGGUAGAACUCACACU	16	24	25885
hsa-miR-222-3p	AGCUACAUCUGGCUACUGGGU	36	21	22187
hsa-miR-125a-5p	UCCCUGAGACCCUUUAACCGUGA	35	24	20960
hsa-miR-16-5p	UAGCAGCACGUAAAUAUUGGCG	29	22	19856
hsa-let-7b-5p	UGAGGUAGUAGGUUGUGUGGUU	28	22	19774
hsa-miR-151a-5p	UCGAGGAGCUCACAGUCUAGU	37	21	19773
hsa-let-7e-5p	UGAGGUAGGAGGUUGUAUAGUU	27	22	19035
hsa-miR-125b-5p	UCCCUGAGACCCUAACUUGUGA	42	22	17965
hsa-let-7i-5p	UGAGGUAGUAGUUUGUGCUGUU	22	22	17802
hsa-let-7g-5p	UGAGGUAGUAGUUUGUACAGUU	43	22	15467
hsa-miR-409-3p	GAAUGUUGCUCGGUGAACCCCU	47	22	14133
hsa-miR-30e-5p	UGUAAACAUCCUUGACUGGAAG	45	22	13889
hsa-miR-181b-5p	AACAUUCAUUGCUGUCGGUGGGU	38	23	12606
hsa-miR-186-5p	CAAAGAAUUCUCCUUUUGGGCU	40	22	12441
hsa-miR-381	UAUACAAGGGCAAGCUCUCUGU	51	22	9851
hsa-miR-423-5p	UGAGGGGCAGAGAGCGAGACUUU	41	23	8893
hsa-miR-30c-5p	UGUAAACAUCCUACACUCUCAGC	66	23	8737
hsa-miR-410	AAUAUAACACAGAUGGCCUGU	50	21	8509
hsa-miR-146b-5p	UGAGAACUGAAUCCAUAAGGCU	19	22	8434
hsa-miR-654-3p	UAUGUCUGCUGACCAUCACCUU	336	22	8392

EP 3 033 418 B1

(continued)

5	CELLS - CTX0E03 07EI				
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
10	hsa-miR-9-5p	UCUUUGGUUAUCUAGCUGUAUGA	58	23	7957
	hsa-miR-28-3p	CACUAGAUUGUGAGCUCCUGGA	56	22	7767
	hsa-miR-148a-3p	UCAGUGCACUACAGAACUUUGU	54	22	6599
	hsa-miR-379-5p	UGGUAGACUAUGGAACGUAGG	18	21	6135
	hsa-let-7d-5p	AGAGGUAGUAGGUUGCAUAGUU	53	22	5972
15	hsa-miR-183-5p	UAUGGCACUGGUAGAAUUCACU	24	22	5477
	hsa-miR-25-3p	CAUUGCACUUGUCUCGGUCUGA	63	22	5303
	hsa-miR-423-3p	AGCUCGGUCUGAGGCCCCUCAGU	57	23	5225
	hsa-miR-889	UUAUAUACGGACAACCAUUGU	64	21	4597
	hsa-miR-221-5p	ACCUGGCAUACAAUGUAGAUUU	39	22	4379
20	hsa-miR-125b-1-3p	ACGGGUUAGGCUCUUGGGAGCU	49	22	4192
	hsa-miR-409-5p	AGGUUACCCGAGCAACUUUGCAU	32	23	3970
	hsa-miR-4492	GGGGCUGGGCGCGCGCC	34	17	3864
	hsa-miR-148b-3p	UCAGUGCAUCACAGAACUUUGU	48	22	3593
	hsa-miR-103a-3p	AGCAGCAUUGUACAGGGCUAUGA	62	23	3518
25	hsa-miR-1271-5p	CUUGGCACCUAGCAAGCACUCA	72	22	3477
	hsa-miR-136-3p	CAUCAUCGUCUCAAUGAGUCU	82	22	3373
	hsa-miR-769-5p	UGAGACCUCUGGGUUCUGAGCU	75	22	2957
	hsa-miR-4532	CCCCGGGGAGCCCGGCG	23	17	2915
	hsa-miR-378a-3p	ACUGGACUUGGAGUCAGAAGG	65	21	2895
30	hsa-miR-99a-5p	AACCCGUAGAUCCGAUCUUGUG	52	22	2767
	hsa-miR-221-3p	AGCUACAUUGUCUGCUGGGUUUC	79	23	2764
	hsa-miR-30e-3p	CUUUCAGUCGGAUGUUUACAGC	71	22	2441
	hsa-miR-26b-5p	UUCAAGUAAUUCAGGAUAGGU	90	21	2432
	hsa-miR-4488	AGGGGGCGGGCUCCGGCG	61	18	2391
35	hsa-miR-27a-3p	UUCACAGUGGCUAAGUUCCGC	46	21	2385
	hsa-miR-23b-3p	AUCACAUUGCCAGGGAUUACC	59	21	2316
	hsa-miR-500a-3p	AUGCACCUGGGCAAGGAUUCUG	44	22	2144
	hsa-miR-941	CACCCGGCUGUGUGCACAUGUGC	60	23	2114
	hsa-miR-23a-3p	AUCACAUUGCCAGGGAUUUCC	55	21	2086
40	hsa-miR-30a-3p	CUUUCAGUCGGAUGUUUGCAGC	77	22	2045
	hsa-miR-30b-5p	UGUAAACAUCCUACACUCAGCU	96	22	1936
	hsa-miR-501-3p	AAUGCACCCGGGCAAGGAUUCU	26	22	1895
	hsa-miR-130b-3p	CAGUGCAAUGAUGAAAGGGCAU	86	22	1862
	hsa-miR-1246	AAUGGAUUUUUGGAGCAGG	21	19	1783

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-140-3p	UACCACAGGGUAGAACCACGG	138	21	1735
hsa-miR-31-5p	AGGCAAGAUGCUGGCAUAGCU	80	21	1705
hsa-miR-493-3p	UGAAGGUCUACUGUGUGCCAGG	83	22	1698
hsa-miR-181c-5p	AACAUUCAACCUGUCGGUGAGU	105	22	1554
hsa-miR-93-5p	CAAAGUGCUGUUCGUGCAGGUAG	116	23	1492
hsa-miR-181a-2-3p	ACCACUGACCGUUGACUGUACC	102	22	1491
hsa-miR-15b-5p	UAGCAGCACAUCAUGGUUUACA	78	22	1465
hsa-miR-7-5p	UGGAAGACUAGUGAUUUUGUUGU	85	23	1460
hsa-miR-192-5p	CUGACCUAUGAAUUGACAGCC	87	21	1453
hsa-miR-425-5p	AAUGACACGAUCACUCCCGUUGA	111	23	1432
hsa-miR-204-5p	UUCCCUUUGUCAUCCUAUGCCU	89	22	1378
hsa-miR-340-5p	UUAUAAAGCAAUGAGACUGAUU	93	22	1360
hsa-miR-190a	UGAU AUGUUUGAUUAUUAGGU	103	22	1305
hsa-miR-34a-5p	UGGCAGUGUCUUAGCUGGUUGU	101	22	1283
hsa-miR-20a-5p	UAAAGUGC UUUAUGUGCAGGUAG	146	23	1257
hsa-miR-29a-3p	UAGCACCAUCUGAAAUCGGUUA	106	22	1206
hsa-miR-361-5p	UUAUCAGAAUCUCCAGGGGUAC	70	22	1173
hsa-miR-671-3p	UCCGGUUCUCAGGGCUCCACC	69	21	1166
hsa-miR-411-5p	UAGUAGACCGUAUAGCGUACG	108	21	1130
hsa-miR-589-5p	UGAGAACCACGUCUGCUCUGAG	73	22	1067
hsa-miR-130a-3p	CAGUGCAAUGU UAAAAGGGCAU	113	22	1020
hsa-miR-320a	AAAAGCUGGGUUGAGAGGGCGA	97	22	994
hsa-miR-149-5p	UCUGGCUCCGUGUCUUCACUCCC	121	23	948
hsa-miR-335-5p	UCAAGAGCAAUAACGAAAAAUGU	128	23	945
hsa-miR-134	UGUGACUGGUUGACCAGAGGGG	94	22	941
hsa-miR-17-5p	CAAAGUGC UUACAGUGCAGGUAG	145	23	939
hsa-miR-493-5p	UUGUACAUGGUAGGCUUUCAUU	88	22	876
hsa-miR-34c-5p	AGGCAGUGUAGUUAGCUGAUUGC	125	23	846
hsa-miR-484	UCAGGCUCAGUCCCCUCCGAU	118	22	835
hsa-miR-181a-3p	ACCAUCGACCGUUGAUUGUACC	100	22	803
hsa-miR-24-3p	UGGCUCAGUUCAGCAGGAACAG	119	22	740
hsa-miR-128	UCACAGUGAACCGGUCUCUUU	109	21	707
hsa-miR-342-3p	UCUCACACAGAAAUCGCACCCGU	81	23	698
hsa-miR-454-3p	UAGUGCAAUAUUGCUUAUAGGGU	169	23	690
hsa-miR-1307-5p	UCGACCGGACCUCGACCGGCU	91	21	616

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-487b	AAUCGUACAGGGUCAUCCACUU	117	22	590
hsa-miR-130b-5p	ACUCUUUCCCGUUGCACUAC	112	21	568
hsa-miR-197-3p	UUCACCACCUUCUCCACCCAGC	122	22	544
hsa-miR-432-5p	UCUUGGAGUAGGUCAUUGGGUGG	95	23	542
hsa-miR-374a-5p	UUAUAAUACAACCUGAUAGUG	74	22	537
hsa-miR-345-5p	GCUGACUCCUAGUCCAGGGCUC	76	22	527
hsa-miR-744-5p	UGCGGGGCUAGGGCUAACAGCA	99	22	515
hsa-miR-376c	AACAUAGAGGAAAUCCACGU	185	21	506
hsa-miR-181d	AACAUUCAUUGUUGUCGGUGGGU	157	23	497
hsa-miR-363-3p	AAUUGCACGGUAUCCAUCUGUA	131	22	493
hsa-miR-539-3p	AUCAUACAAGGACAAUUCUUU	150	22	493
hsa-miR-758	UUUGUGACCUGGUCCACUAACC	141	22	477
hsa-miR-323a-3p	CACAUUACACGGUCGACCUCU	158	21	443
hsa-miR-107	AGCAGCAUUGUACAGGGCUAUCA	254	23	431
hsa-miR-720	UCUCGCUGGGGCCUCCA	84	17	427
hsa-miR-654-5p	UGGUGGGCCGCAGAACAUUGUC	115	22	409
hsa-miR-370	GCCUGCUGGGGUGGAACCUGGU	126	22	406
hsa-miR-421	AUCAACAGACAUUAAUUGGGCGC	151	23	399
hsa-miR-30d-3p	CUUUCAGUCAGAUUUUGCUGC	114	22	358
hsa-miR-148b-5p	AAGUUCUGUUAUACACUCAGGC	127	22	354
hsa-miR-1301	UUGCAGCUGCCUGGGAGUGACUUC	181	24	346
hsa-miR-374b-5p	AUAUAAUACAACCUGCUAAGUG	143	22	339
hsa-miR-125b-2-3p	UCACAAGUCAGGCUCUUGGGAC	68	22	333
hsa-miR-28-5p	AAGGAGCUCACAGUCUAUUGAG	152	22	332
hsa-miR-495	AAACAAACAUGGUGCACUUCUU	241	22	321
hsa-miR-15a-5p	UAGCAGCACAUAAUGGUUUGUG	223	22	320
hsa-miR-100-3p	CAAGCUUGUAUCUAUAGGUAUG	98	22	314
hsa-miR-193b-3p	AACUGGCCCUCAAAGUCCCGCU	148	22	305
hsa-miR-330-5p	UCUCUGGGCCUGUGUCUUAGGC	161	22	303
hsa-miR-376a-3p	AUCAUAGAGGAAAAUCCACGU	237	21	298
hsa-miR-135b-5p	UAUGGCUUUUCAUUCCUAUGUGA	137	23	289
hsa-miR-301a-3p	CAGUGCAAUAGUAUUGUCAAGC	107	23	280
hsa-miR-218-5p	UUGUGCUUGAUCUAACCAUGU	206	21	276
hsa-miR-143-3p	UGAGAUGAAGCACUGUAGCUC	176	21	256
hsa-miR-27b-5p	AGAGCUUAGCUGAUUGGUGAAC	201	22	255

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-369-3p	AAUAAUACAUGGUUGAUCUUU	196	21	255
hsa-miR-877-5p	GUAGAGGAGAUGGCGCAGGG	133	20	249
hsa-miR-19b-3p	UGUGCAAAUCCAUGCAAAACUGA	163	23	246
hsa-miR-424-5p	CAGCAGCAAUUCAUGUUUUGAA	186	22	245
hsa-miR-660-5p	UACCCAUUGCAUAUCGGAGUUG	187	22	244
hsa-miR-532-5p	CAUGCCUUGAGUGUAGGACCGU	178	22	238
hsa-miR-299-3p	UAUGUGGGAUGGUAACCGCUU	182	22	235
hsa-miR-431-3p	CAGGUCGUCUUGCAGGGCUUCU	168	22	231
hsa-miR-374a-3p	CUUAUCAGAUUGUAUUGUAAUU	173	22	220
hsa-miR-148a-5p	AAAGUUCUGAGACACUCCGACU	144	22	214
hsa-miR-4516	GGGAGAAGGGUCGGGGC	110	17	207
hsa-miR-92b-5p	AGGGACGGGACGCGGUGCAGUG	208	22	206
hsa-miR-16-2-3p	CCAAUAUUACUGUGCUGCUUUA	316	22	202
hsa-miR-101-3p	UACAGUACUGUGAUAACUGAA	142	21	201
hsa-let-7a-3p	CUAUACAAUCUACUGUCUUUC	222	21	199
hsa-miR-4485	UAACGGCCGCGGUACCCUAA	67	20	195
hsa-miR-455-3p	GCAGUCCAUGGGCAUAUACAC	140	21	192
hsa-miR-185-5p	UGGAGAGAAAGGCAGUUCCUGA	224	22	188
hsa-miR-1185-1-3p	AUAUACAGGGGGAGACUCUUAU	209	22	187
hsa-miR-1197	UAGGACACAUGGUCUACUUCU	244	21	185
hsa-miR-106b-3p	CCGCACUGUGGGUACUUGCUGC	159	22	178
hsa-miR-24-2-5p	UGCCUACUGAGCUGAAACACAG	156	22	178
hsa-miR-4677-3p	UCUGUGAGACCAAAGAACUACU	120	22	177
hsa-miR-380-3p	UAUGUAAUAUGGUCCACAUCUU	445	22	174
hsa-miR-548k	AAAAGUACUUGCGGAUUUUGCU	198	22	171
hsa-miR-1307-3p	ACUCGGCGUGGGCUCGGUCGUG	124	22	169
hsa-miR-485-3p	GUCAUACACGGCUCUCCUCUCU	153	22	168
hsa-miR-494	UGAAACAUAACACGGGAAACCUC	240	22	165
hsa-miR-17-3p	ACUGCAGUGAAGGCACUUGUAG	184	22	163
hsa-miR-561-5p	AUCAAGGAUCUUAACUUUGCC	179	22	160
hsa-miR-27a-5p	AGGGCUUAGCUGCUUGUGAGCA	130	22	158
hsa-miR-874	CUGCCCUGGCCCGAGGGACCGA	147	22	151
hsa-miR-9-3p	AUAAAGCUAGAUAAACCGAAAGU	183	22	151
hsa-miR-96-5p	UUUGGCACUAGCACAUUUUUGCU	123	23	151
hsa-miR-656	AAUAUUUAUCAGUCAACCUCU	221	21	147

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-379-3p	UAUGUAACAUGGUCCACUAACU	446	22	145
hsa-miR-382-5p	GAAGUUGUUCGUGGUGGAUUCG	238	22	144
hsa-miR-541-3p	UGGUGGGGCACAGAAUCUGGACU	136	22	141
hsa-miR-337-3p	CUCCUAUAUGAUGCCUUUCUUC	215	22	139
hsa-miR-15b-3p	CGAAUCAUUAUUUGCUGCUCUA	447	22	137
hsa-miR-20b-5p	CAAAGUGCUCAUAGUGCAGGUAG	317	23	136
hsa-miR-329	AACACACCUGGUUAACCUCUUU	214	22	136
hsa-miR-3676-5p	AGGAGAUCCUGGGUU	280	15	134
hsa-miR-543	AAACAUUCGCGGUGCACUUCUU	193	22	134
hsa-miR-365b-3p	UAAUGCCCCUAAAAUCCUUAU	279	22	133
hsa-miR-125a-3p	ACAGGUGAGGUUCUUGGGAGCC	160	22	131
hsa-miR-3065-5p	UCAACAAAUCACUGAUGCUGGA	226	23	130
hsa-miR-1296	UUAGGGCCCUGGCUCCAUCUCC	271	22	126
hsa-miR-935	CCAGUUACCGCUUCCGCUACCGC	311	23	118
hsa-miR-132-3p	UACAGUCUACAGCCAUGGUCG	104	22	116
hsa-miR-4284	GGGCUCACAUCACCCCAU	191	18	116
hsa-miR-487a	AAUCAUACAGGGACAUCAGUU	203	22	113
hsa-miR-574-5p	UGAGUGUGUGUGUGUGAGUGUGU	334	23	113
hsa-miR-301b	CAGUGCAAUGAUUUGUCAAGC	164	23	111
hsa-miR-548o-3p	CCAAAACUGCAGUUACUUUUGC	268	22	105
hsa-miR-18a-5p	UAAGGUGCAUCUAGUGCAGAUAG	256	23	104
hsa-miR-485-5p	AGAGGCUGGCCGUGAUGAAUUC	165	22	104
hsa-miR-548ah-5p	AAAAGUGAUUGCAGUGUUUG	235	20	103
hsa-miR-361-3p	UCCCCAGGUGUGAUUCUGAUUU	250	23	101
hsa-miR-433	AUCAUGAUGGGCUCCUCGGUGU	174	22	101
hsa-miR-337-5p	GAACGGCUUCAUACAGGAGUU	277	21	100
hsa-miR-1276	UAAAGAGCCCUGUGGAGACA	194	20	99
hsa-miR-30c-1-3p	CUGGGAGAGGGUUGUUUACUCC	213	22	99
hsa-miR-31-3p	UGCUAUGCCAACAUAUUGCCA	172	22	96
hsa-miR-424-3p	CAAAACGUGAGGCGCUGCUAU	298	21	96
hsa-miR-550a-5p	AGUGCCUGAGGGAGUAAGAGCCC	134	23	95
hsa-miR-4454	GGAUCCGAGUCACGGCACCA	299	20	94
hsa-miR-541-5p	AAAGGAUUCUGCUGUCGGUCCACU	432	25	92
hsa-miR-106b-5p	UAAAGUGCUGACAGUGCAGAU	170	21	89
hsa-miR-153	UUGCAUAGUCACAAAAGUGAUC	188	22	88

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-135b-3p	AUGUAGGGCUAAAAGCCAUGGG	205	22	87
hsa-miR-574-3p	CACGCUCAUGCACACACCCACA	253	22	87
hsa-miR-1226-3p	UCACCAGCCCUGUGUUCCCUAG	199	22	85
hsa-miR-576-5p	AUUCUAAUUUCUCCACGUCUUU	306	22	84
hsa-miR-127-5p	CUGAAGCUCAGAGGGCUCUGAU	255	22	83
hsa-miR-155-5p	UUAAUUGCUAUUCGUGAUAGGGGU	448	23	83
hsa-miR-3176	ACUGGCCUGGGACUACCGG	227	19	83
hsa-miR-382-3p	AAUCAUUCACGGACAACACUU	260	21	83
hsa-miR-1275	GUGGGGGAGAGGCUGUC	162	17	82
hsa-miR-671-5p	AGGAAGCCCUGGAGGGGCUGGAG	288	23	82
hsa-miR-23a-5p	GGGGUUCCUGGGGAUGGGAUUU	212	22	81
hsa-miR-25-5p	AGGCGGAGACUUGGGCAAUUG	225	21	80
hsa-miR-641	AAAGACAUAGGAUAGAGUCACCUC	285	24	80
hsa-miR-19a-3p	UGUGCAAAUCUAUGCAAAACUGA	177	23	79
hsa-miR-377-3p	AUCACACAAAGGCAACUUUUGU	449	22	78
hsa-miR-454-5p	ACCCUAUCAUAUUGUCUCUGC	265	22	78
hsa-miR-496	UGAGUAUUACAUGGCCAAUCUC	267	22	78
hsa-miR-29b-3p	UAGCACCAUUUGAAAUCAGUGUU	166	23	77
hsa-miR-26a-2-3p	CCUAUUCUUGAUUACUUGUUUC	257	22	76
hsa-miR-1260b	AUCCCACCACUGCCACCAU	373	19	74
hsa-miR-2467-5p	UGAGGCUCUGUUAGCCUUGGCUC	154	23	74
hsa-miR-377-5p	AGAGGUUGCCCUUGGUGAAUUC	202	22	74
hsa-miR-330-3p	GCAAAGCACACGGCCUGCAGAGA	195	23	73
hsa-miR-1180	UUUCCGGCUCGCGUGGGUGUGU	313	22	71
hsa-miR-99b-3p	CAAGCUCGUGUCUGUGGGUCCG	243	22	71
hsa-miR-299-5p	UGGUUUACCGUCCCACAUACAU	319	22	69
hsa-miR-374b-3p	CUUAGCAGGUUGUAUUAUCAUU	229	22	69
hsa-miR-4746-5p	CCGGUCCCAGGAGAACCUGCAGA	266	23	69
hsa-miR-331-3p	GCCCCUGGGCCUAUCCUAGAA	450	21	68
hsa-miR-340-3p	UCCGUCUCAGUUACUUUAUAGC	248	22	68
hsa-miR-92a-1-5p	AGGUUGGGAUCGGUUGCAAUGCU	204	23	68
hsa-miR-542-3p	UGUGACAGAUUGAUAAACUGAAA	331	22	66
hsa-miR-431-5p	UGUCUUGCAGGCCGUAUGCA	132	21	65
hsa-miR-1254	AGCCUGGAAGCUGGAGCCUGCAGU	270	24	61
hsa-miR-3158-3p	AAGGGCUUCCUCUCUGCAGGAC	167	22	61

(continued)

5	CELLS - CTX0E03 07EI				
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
10	hsa-miR-362-5p	AAUCCUUGGAACCUAGGUGUGAGU	139	24	61
	hsa-miR-30c-2-3p	CUGGGAGAAGGCUGUUUACUCU	321	22	59
	hsa-miR-4461	GAUUGAGACUAGUAGGGCUAGGC	129	23	59
	hsa-miR-3200-3p	CACCUUGCGCUACUCAGGUCUG	247	22	57
	hsa-miR-215	AUGACCUAUGAAUUGACAGAC	451	21	56
15	hsa-miR-1185-5p	AGAGGAUACCCUUUGUAUGUU	368	21	55
	hsa-miR-328	CUGGCCUCUCUGCCCUUCCGU	297	22	55
	hsa-miR-655	AUAAUACAUGGUUAACCUCUUU	286	22	55
	hsa-miR-181b-3p	CUCACUGAACAAUGAAUGCAA	245	21	54
	hsa-miR-376b	AUCAUAGAGGAAAAUCCAUGUU	452	22	54
20	hsa-miR-486-3p	CGGGGCAGCUCAGUACAGGAU	453	21	54
	hsa-miR-760	CGGCUCUGGGUCUGUGGGGA	289	20	54
	hsa-miR-3909	UGUCCUCUAGGGCCUGCAGUCU	412	22	53
	hsa-miR-4508	GCGGGGCUGGGCGCGCG	135	17	53
	hsa-miR-4521	GCUAAGGAAGUCCUGUGCUCAG	233	22	53
25	hsa-let-7e-3p	CUAUACGGCCUCCUAGCUUCC	290	22	52
	hsa-miR-455-5p	UAUGUGCCUUUGGACUACAUCG	192	22	52
	hsa-miR-93-3p	ACUGCUGAGCUAGCACUCCCCG	454	22	51
	hsa-miR-151b	UCGAGGAGCUCACAGUCU	455	18	49
	hsa-miR-887	GUGAACGGGCGCCAUCCCGAGG	456	22	49
30	hsa-miR-152	UCAGUGCAUGACAGAACUUGG	344	21	48
	hsa-miR-324-3p	ACUGCCCCAGGUGCUGCUGG	276	20	48
	hsa-miR-1266	CCUCAGGGCUGUAGAACAGGGCU	457	23	47
	hsa-miR-302b-3p	UAAGUGCUUCCAUGUUUUAGUAG	458	23	47
	hsa-miR-548e	AAAAACUGAGACUACUUUUGCA	459	22	47
35	hsa-miR-502-3p	AAUGCACCUGGGCAAGGAUUCA	281	22	46
	hsa-miR-302d-3p	UAAGUGCUUCCAUGUUUGAGUGU	460	23	45
	hsa-miR-3943	UAGCCCCCAGGCUUCACUUGGCG	207	23	45
	hsa-miR-1286	UGCAGGACCAAGAUGAGCCCU	293	21	44
	hsa-miR-3605-5p	UGAGGAUGGAUAGCAAGGAAGCC	189	23	44
40	hsa-miR-505-3p	CGUCAACACUUGCUGGUUUCCU	282	22	44
	hsa-miR-3615	UCUCUCGGCUCCUCGCGGCUC	323	21	43
	hsa-miR-4435	AUGGCCAGAGCUCACACAGAGG	230	22	43
	hsa-miR-598	UACGUCAUCGUUGUCAUCGUCA	461	22	43
	hsa-miR-126-5p	CAUUAUUACUUUUGGUACGCG	462	21	42

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-4671-3p	UUAGUGCAUAGUCUUUGGUCU	301	21	41
hsa-miR-652-3p	AAUGGCGCCACUAGGGUUGUG	442	21	41
hsa-miR-3687	CCCGGACAGGCGUUCGUGCGACGU	190	24	40
hsa-miR-4286	ACCCACUCCUGGUACC	328	17	40
hsa-miR-590-3p	UAAUUUUAUGUAUAAGCUAGU	463	21	40
hsa-miR-1285-3p	UCUGGGCAACAAAGUGAGACCU	464	22	39
hsa-miR-2355-5p	AUCCCCAGAUACAAUGGACAA	593	21	38
hsa-miR-550a-3p	UGUCUUACUCCCUCAGGCACAU	283	22	38
hsa-let-7d-3p	CUAUACGACCUGCUGCCUUUCU	92	22	37
hsa-miR-136-5p	ACUCCAUUUGUUUUGAUGAUGGA	272	23	37
hsa-miR-1468	CUCCGUUUGCCUGUUUCGCUG	296	21	37
hsa-miR-3609	CAAAGUGAUGAGUAAUACUGGCUG	216	24	37
hsa-miR-548b-5p	AAAAGUAAUUGUGGUUUUGGCC	304	22	37
hsa-miR-664-3p	UAUUCAUUUAUCCCCAGCCUACA	287	23	37
hsa-miR-99a-3p	CAAGCUCGCUUCUAUGGGUCUG	367	22	37
hsa-miR-532-3p	CCUCCCACACCCAAGGCUUGCA	252	22	36
hsa-miR-10b-5p	UACCCUGUAGAACCGAAUUUGUG	465	23	33
hsa-miR-369-5p	AGAUCGACCGUGUUUAUUCGC	357	22	33
hsa-miR-3161	CUGAUAAGAACAGAGGCCCAGAU	466	23	32
hsa-miR-3940-3p	CAGCCCGGAUCCCAGCCCACUU	239	22	32
hsa-miR-663b	GGUGGCCCGGCCGUGCCUGAGG	180	22	32
hsa-miR-219-2-3p	AGAAUUGUGGCUGGACAUCUGU	467	22	31
hsa-miR-2277-5p	AGCGCGGGCUGAGCGCUGCCAGUC	735	24	31
hsa-miR-4448	GGCUCCUUGGUCUAGGGGUA	231	20	31
hsa-miR-339-5p	UCCCUGUCCUCCAGGAGCUCACG	402	23	30
hsa-miR-3613-5p	UGUUGUACUUUUUUUUUGUUC	469	22	30
hsa-miR-4775	UUAAUUUUUUGUUUCGGUCACU	302	22	30
hsa-miR-212-5p	ACCUUGGCUCUAGACUGCUUACU	246	23	29
hsa-miR-324-5p	CGCAUCCCCUAGGGCAUUGGUGU	354	23	27
hsa-miR-4326	UGUUCCUCUGUCUCCCAGAC	417	20	27
hsa-miR-582-3p	UACUGGUUGAACAACUGAACC	470	22	27
hsa-miR-34a-3p	CAAUCAGCAAGUAUACUGCCCU	403	22	26
hsa-miR-106a-5p	AAAAGUGCUUACAGUGCAGGUAG	471	23	25
hsa-miR-4745-5p	UGAGUGGGGCUCCCGGGACGGCG	219	23	25
hsa-miR-769-3p	CUGGGAUCUCCGGGGUCUUGGUU	337	23	25

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-1268a	CGGGCGUGGUGGUGGGGG	291	18	24
hsa-miR-154-3p	AAUCAUACACGGUUGACCUAUU	472	22	24
hsa-miR-188-3p	CUCCCACAUGCAGGGUUUGCA	200	21	24
hsa-miR-29c-3p	UAGCACCAUUUGAAAUCGGUUA	473	22	24
hsa-miR-539-5p	GGAGAAAUUAUCCUUGGUGUGU	234	22	24
hsa-miR-766-3p	ACUCCAGCCCCACAGCCUCAGC	310	22	24
hsa-miR-30b-3p	CUGGGAGGUGGAUGUUUACUUC	320	22	23
hsa-miR-3177-3p	UGCACGGCACUGGGGACACGU	275	21	23
hsa-miR-191-3p	GCUGCGCUUGGAUUUCGUCCCC	474	22	22
hsa-miR-296-3p	GAGGGUUGGGUGGAGGCUCUCC	274	22	22
hsa-miR-296-5p	AGGGCCCCCCCUCAAUCCUGU	258	21	22
hsa-miR-339-3p	UGAGCGCCUCGACGACAGAGCCG	228	23	22
hsa-miR-501-5p	AAUCCUUUGUCCCUGGGUGAGA	430	22	22
hsa-miR-200b-3p	UAAUACUGCCUGGUAAUGAUGA	475	22	21
hsa-miR-212-3p	UAACAGUCUCCAGUCACGGCC	348	21	21
hsa-miR-26b-3p	CCUGUUCUCCAUAUACUUGGCUC	391	22	21
hsa-miR-665	ACCAGGAGGCUGAGGCCCCU	309	20	21
hsa-miR-668	UGUCACUCGGCUCGGCCCACUAC	476	23	21
hsa-miR-146a-5p	UGAGAACUGAAUCCAUGGGUU	477	22	20
hsa-miR-1973	ACCGUGCAAAGGUAGCAUA	171	19	20
hsa-miR-210	CUGUGCGUGUGACAGCGGCUGA	478	22	20
hsa-miR-3607-5p	GCAUGUGAUGAAGCAAUUCAGU	249	22	20
hsa-miR-378a-5p	CUCCUGACUCCAGGUCCUGUGU	217	22	20
hsa-miR-4449	CGUCCCGGGGCUGCGCGAGGCA	155	22	20
hsa-miR-138-5p	AGCUGGUGUUGUGAAUCAGGCCG	379	23	19
hsa-miR-146b-3p	UGCCCUGUGGACUCAGUUCUGG	381	22	18
hsa-miR-3065-3p	UCAGCACCAGGAUAUUGUUGGAG	350	23	18
hsa-miR-4417	GGUGGGCUUCCCGGAGGG	175	18	18
hsa-miR-497-5p	CAGCAGCACACUGUGGUUUUGU	479	21	18
hsa-miR-500a-5p	UAAUCCUUGCUACCUGGGUGAGA	303	23	18
hsa-miR-625-3p	GACUAUAGAACUUUCCCCUCA	307	22	18
hsa-miR-628-3p	UCUAGUAAGAGUGGCAGUCGA	335	21	18
hsa-miR-1343	CUCCUGGGGCCCCGCACUCUCGC	378	22	17
hsa-miR-3648	AGCCGCGGGGAUCGCCGAGGG	259	21	17
hsa-miR-432-3p	CUGGAUGGCUCCUCCAUGUCU	262	21	17

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-4482-3p	UUUCUAAUUUCUCAGUGGGGCUC	361	22	17
hsa-miR-542-5p	UCGGGGAUCAUCAUGUCACGAGA	433	23	17
hsa-miR-551b-3p	GCGACCCAUACUUGGUUUCAG	434	21	17
hsa-miR-7-1-3p	CAACAAAUACACAGUCUGCCAU	480	22	17
hsa-miR-219-1-3p	AGAGUUGAGUCUGGACGUCCCG	390	22	16
hsa-miR-3656	GGCGGGUGCGGGGGUGG	251	17	16
hsa-miR-3661	UGACCUGGGACUCGGACAGCUG	481	22	16
hsa-miR-411-3p	UAUGUAACACGGUCCACUAACC	482	22	16
hsa-miR-5096	GUUUCACCAUGUUGGUCAGGC	220	21	16
hsa-miR-577	UAGAUAAAAUAUUGGUACCUG	436	21	16
hsa-let-7i-3p	CUGCGCAAGCUACUGCCUUGCU	483	22	15
hsa-miR-132-5p	ACCGUGGCUUUCGAUUGUUACU	315	22	15
hsa-miR-140-5p	CAGUGGUUUUACCCUAUGGUAG	380	22	15
hsa-miR-195-5p	UAGCAGCACAGAAUAUUGGC	346	21	15
hsa-miR-3187-3p	UUGGCCAUGGGGCUGCGCGG	322	20	15
hsa-miR-342-5p	AGGGGUGCUAUCUGUGAUUGA	278	21	15
hsa-miR-34b-3p	CAAUCACUAACUCCACUGCCAU	404	22	15
hsa-miR-4661-5p	AACUAGCUCUGUGGAUCCUGAC	484	22	15
hsa-miR-584-5p	UUAUGGUUUUGCCUGGGACUGAG	485	22	15
hsa-miR-744-3p	CUGUUGCCACUAACCUCACCCU	486	22	15
hsa-miR-770-5p	UCCAGUACCACGUGUCAGGGCCA	487	23	15
hsa-miR-3677-3p	CUCGUGGGCUCUGGCCACGGCC	356	22	14
hsa-miR-425-3p	AUCGGGAAUGUCGUGUCCGCCC	358	22	14
hsa-miR-548ah-3p	CAAAAACUGCAGUUACUUUUGC	149	22	14
hsa-miR-5699	UCCUGUCUUUCCUUGUUGGAGC	488	22	14
hsa-miR-582-5p	UUACAGUUGUUAACCAGUUACU	489	23	14
hsa-miR-1185-2-3p	AUAUACAGGGGGAGACUCUCAU	314	22	13
hsa-miR-1249	ACGCCCUUCSCCCCCCUUCUUA	490	22	13
hsa-miR-1255a	AGGAUGAGCAAAGAAAGUAGAUU	341	23	13
hsa-miR-1910	CCAGUCCUGUGCCUGCCGCCU	236	21	13
hsa-miR-301a-5p	GCUCUGACUUUAUUGCACUACU	491	22	13
hsa-miR-5001-3p	UUCUGCCUCUGUCCAGGUCCUU	492	22	13
hsa-miR-5094	AAUCAGUGAAUGCCUUGAACCU	493	22	13
hsa-miR-628-5p	AUGCUGACAUUUUACUAGAGG	440	22	13
hsa-miR-629-5p	UGGGUUUACGUUGGGAGAACU	441	21	13

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-937	AUCCGCGCUCUGACUCUCUGCC	312	22	13
hsa-miR-940	AAGGCAGGGCCCCCGCUCCCC	366	21	13
hsa-miR-1248	ACCUUCUUGUAUAAGCACUGUGCUGAAA	269	27	12
hsa-miR-194-5p	UGUAACAGCAACUCCAUGUGGA	345	22	12
hsa-miR-199b-3p	ACAGUAGUCUGCACAUUGGUUA	494	22	12
hsa-miR-22-5p	AGUUCUUCAGUGGCAAGCUUUA	495	22	12
hsa-miR-3605-3p	CCUCCGUGUUACCUGUCCUCUAG	496	23	12
hsa-miR-3654	GACUGGACAAGCUGAGGAA	325	19	12
hsa-miR-504	AGACCCUGGUCUGCACUCUAUC	497	22	12
hsa-miR-1291	UGGCCCUGACUGAAGACCAGCAGU	294	24	11
hsa-miR-1299	UUCUGGAAUUCUGUGUGAGGGA	498	22	11
hsa-miR-188-5p	CAUCCCUUGCAUGGUGGAGGG	499	21	11
hsa-miR-222-5p	CUCAGUAGCCAGUGUAGAUCU	349	22	11
hsa-miR-331-5p	CUAGGUAUGGUCCCAGGGAUCC	500	22	11
hsa-miR-3939	UACGCGCAGACCACAGGAUGUC	261	22	11
hsa-miR-154-5p	UAGGUUAUCCGUGUUGCCUUCG	501	22	10
hsa-miR-18a-3p	ACUGCCCUAAGUGCUCUCCUUCUGG	502	23	10
hsa-miR-1908	CGGCGGGGACGGCGAUUGGUC	383	21	10
hsa-miR-200c-3p	UAAUACUGCCGGGUAAUGAUGGA	347	23	10
hsa-miR-2116-3p	CCUCCCAUGCCAAGAACUCCC	318	21	10
hsa-miR-302a-3p	UAAGUGCUUCCAUGUUUUGGUGA	503	23	10
hsa-miR-3174	UAGUGAGUUAGAGAUGCAGAGCC	353	23	10
hsa-miR-326	CCUCUGGGCCCUUCCUCCAG	504	20	10
hsa-let-7g-3p	CUGUACAGGCCACUGCCUUGC	505	21	9
hsa-miR-141-3p	UAACACUGUCUGGUAAAAGAUGG	295	22	9
hsa-miR-24-1-5p	UGCCUACUGAGCUGAUUAUCAGU	506	22	9
hsa-miR-3115	AUAUGGGUUUACUAGUUGGU	351	20	9
hsa-miR-3180-3p	UGGGGCGGAGCUUCCGGAGGCC	400	22	9
hsa-miR-33a-5p	GUGCAUUGUAGUUGCAUUGCA	355	21	9
hsa-miR-34c-3p	AAUCACUAACCACACGGCCAGG	405	22	9
hsa-miR-3929	GAGGCUGAUGUGAGUAGACCACU	218	23	9
hsa-miR-4517	AAAUAUGAUGAAACUCACAGCUGAG	507	25	9
hsa-miR-576-3p	AAGAUGUGGAAAAAUUGGAAUC	508	22	9
hsa-miR-1229	CUCUCACCACUGCCCUCCACAG	509	23	8
hsa-miR-1289	UGGAGUCCAGGAAUCUGCAUUUU	343	23	8

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-1915-5p	ACCUUGCCUUGCUGCCCGGGCC	385	22	8
hsa-miR-23b-5p	UGGGUUCCUGGCAUGCUGAUUU	510	22	8
hsa-miR-302a-5p	ACUUAACGUGGAUGUACUUGCU	511	23	8
hsa-miR-3938	AAUUCCCUUGUAGAUAAACCCGG	512	22	8
hsa-miR-4466	GGGUGCGGGCCGGCGGGG	264	18	8
hsa-miR-4786-5p	UGAGACCAGGACUGGAUGCACC	197	22	8
hsa-miR-589-3p	UCAGAACAAAUGCCGGUUCCAGA	438	24	8
hsa-miR-616-5p	ACUCAAAACCCUUCAGUGACUU	439	22	8
hsa-miR-943	CUGACUGUUGCCGUCCUCCAG	338	21	8
hsa-miR-1237	UCCUUCUGCUCCGUCCCCCAG	370	21	7
hsa-miR-1915-3p	CCCCAGGGCGACGCGGCGGG	384	20	7
hsa-miR-3620	UCACCCUGCAUCCCGCACCCAG	324	22	7
hsa-miR-3691-5p	AGUGGAUGAUGGAGACUCGGUAC	409	23	7
hsa-miR-4426	GAAGAUGGACGUACUUU	359	17	7
hsa-let-7a-2-3p	CUGUACAGCCUCCUAGCUUUC	513	22	6
hsa-miR-10a-3p	CAAAUUCGUAUCUAGGGGAAUA	514	22	6
hsa-miR-1287	UGCUGGAUCAGUGGUUCGAGUC	515	22	6
hsa-miR-145-5p	GUCCAGUUUUCCAGGAAUCCCU	516	23	6
hsa-miR-29b-1-5p	GCUGGUUUAUAUGGUGGUUAGA	517	24	6
hsa-miR-3128	UCUGGCAAGUAAAAACUCUCAU	518	23	6
hsa-miR-33b-5p	GUGCAUUGCUGUUGCAUUGC	519	20	6
hsa-miR-3681-5p	UAGUGGAUGAUGCACUCUGUGC	327	22	6
hsa-miR-3685	UUUCCUACCCUACCUGAAGACU	520	22	6
hsa-miR-3918	ACAGGGCCGCAGAUGGAGACU	521	21	6
hsa-miR-551b-5p	GAAAUCAAGCGUGGGUGAGACC	522	22	6
hsa-miR-1273f	GGAGAUGGAGGUUGCAGUG	292	19	5
hsa-miR-1273g-3p	ACCACUGCACUCCAGCCUGAG	210	21	5
hsa-miR-1304-5p	UUUGAGGCUACAGUGAGAUGUG	523	22	5
hsa-miR-1538	CGGCCCCGGCUGCUGCUGUCCU	524	23	5
hsa-miR-181c-3p	AACCAUCGACCGUUGAGUGGAC	525	22	5
hsa-miR-193a-5p	UGGGUCUUUGCGGGCGAGAUGA	526	22	5
hsa-miR-208b	AUAAGACGAACAAAGGUUUGU	388	22	5
hsa-miR-219-5p	UGAUUGUCCAAACGCAAUUCU	527	21	5
hsa-miR-3159	UAGGAUUACAAGUGUCGGCCAC	528	22	5
hsa-miR-3173-5p	UGCCCUGCCUGUUUUCUCCUUU	529	22	5

EP 3 033 418 B1

(continued)

	CELLS - CTX0E03 07EI				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-3175	CGGGGAGAGAACGCAGUGACGU	530	22	5
	hsa-miR-3200-5p	AAUCUGAGAAGGCGCACAAGGU	531	22	5
10	hsa-miR-3662	GAAAAUGAUGAGUAGUGACUGAUG	326	24	5
	hsa-miR-3928	GGAGGAACCUUGGAGCUUCGGC	413	22	5
	hsa-miR-4709-3p	UUGAAGAGGAGGUGCUCUGUAGC	532	23	5
	hsa-miR-4787-3p	GAUGCGCCGCCACUGCCCCGCGC	533	24	5
	hsa-miR-499a-5p	UUAAGACUUGCAGUGAUGUUU	534	21	5
	hsa-miR-545-3p	UCAGCAAACAUUUUUGUGUGC	242	22	5
	hsa-miR-548u	CAAAGACUGCAAUUACUUUUGCG	535	23	5
20	hsa-miR-659-5p	AGGACCUUCCCUGAACCAAGGA	364	22	5
	hsa-miR-1257	AGUGAAUGAUGGGUUCUGACC	372	21	4
	hsa-miR-1292	UGGGAACGGGUUCCGGCAGACGCUG	536	25	4
25	hsa-miR-1914-5p	CCCUGUGCCCGGCCACUUCUG	537	22	4
	hsa-miR-195-3p	CCAAUAUUGGCUGUGCUGCUCC	538	22	4
	hsa-miR-2110	UUGGGGAAACGGCCGCUGAGUG	389	22	4
	hsa-miR-302c-5p	UUUAACAUGGGGGUACCUGCUG	539	22	4
30	hsa-miR-3126-3p	CAUCUGGCAUCCGUCACACAGA	394	22	4
	hsa-miR-3126-5p	UGAGGGACAGAUGCCAGAAGCA	352	22	4
	hsa-miR-3150a-5p	CAACCUCGACGAUCCUCAGC	540	22	4
35	hsa-miR-3157-3p	CUGCCCUAGUCUAGCUGAAGCU	399	22	4
	hsa-miR-323b-3p	CCCAAUACACGGUCGACCUCUU	541	22	4
	hsa-miR-335-3p	UUUUUCAUUAUUGCUCUGACC	542	22	4
	hsa-miR-3607-3p	ACUGUAAACGCUUUCUGAUG	543	20	4
40	hsa-miR-3653	CUAAGAAGUUGACUGAAG	544	18	4
	hsa-miR-3663-3p	UGAGCACCACACAGGCCGGGCGC	545	23	4
	hsa-miR-376a-5p	GUAGAUUCUCCUUCUAUGAGUA	410	22	4
45	hsa-miR-4423-3p	AUAGGCACCAAAAAGCAACAA	662	21	4
	hsa-miR-4423-5p	AGUUGCCUUUUUGUUCCCAUGC	263	22	4
	hsa-miR-4463	GAGACUGGGGUGGGGCC	300	17	4
	hsa-miR-449a	UGGCAGUGUAUUGUUAGCUGGU	547	22	4
50	hsa-miR-4511	GAAGAACUGUUGCAUUUGCCCU	548	22	4
	hsa-miR-4640-3p	CACCCCUUGUUUCCUGGCCAC	329	22	4
	hsa-miR-4800-3p	CAUCCGUCCGUCUGUCCAC	549	19	4
55	hsa-miR-505-5p	GGGAGCCAGGAAGUAUUGAUGU	550	22	4
	hsa-miR-548a-3p	CAAAACUGGCAAUUACUUUUGC	551	22	4

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-570-3p	CGAAAACAGCAAUUACCUUUGC	333	22	4
hsa-miR-663a	AGGCGGGGCGCCGCGGGACCGC	365	22	4
hsa-miR-877-3p	UCCUCUUCUCCCUCCUCCAG	552	21	4
hsa-miR-103a-2-5p	AGCUUCUUUACAGUGCUGCCUUG	553	23	3
hsa-miR-1268b	CGGGCGUGGUGGUGGGGGUG	554	20	3
hsa-miR-1270	CUGGAGAUUAUGGAAGAGCUGUGU	555	23	3
hsa-miR-1293	UGGGUGGUCUGGAGAUUUGUGC	556	22	3
hsa-miR-1322	GAUGAUGCUGCUGAUGCUG	557	19	3
hsa-miR-150-5p	UCUCCCAACCCUUGUACCAGUG	558	22	3
hsa-miR-190b	UGAUUAUGUUUGAUUUGGGUU	559	21	3
hsa-miR-193a-3p	AACUGGCCUACAAAGUCCAGU	386	22	3
hsa-miR-193b-5p	CGGGGUUUUGAGGGCGAGAUGA	560	22	3
hsa-miR-199a-5p	CCCAGUGUUCAGACUACCUGUUC	273	23	3
hsa-miR-20a-3p	ACUGCAUUAUGAGCACUAAAAG	561	22	3
hsa-miR-216a	UAAUCUCAGCUGGCAACUGUGA	562	22	3
hsa-miR-2682-5p	CAGGCAGUGACUGUUCAGACGUC	563	23	3
hsa-miR-2964a-5p	AGAUGUCCAGCCACAAUUCUCG	564	22	3
hsa-miR-3177-5p	UGUGUACACACGUGCCAGGCGCU	565	23	3
hsa-miR-320c	AAAAGCUGGGUUGAGAGGGU	566	20	3
hsa-miR-323a-5p	AGGUGGUCCGUGGCGGUUCGC	567	22	3
hsa-miR-3622a-5p	CAGGCACGGGAGCUCAGGUGAG	568	22	3
hsa-miR-3912	UAACGCAUAAUAUGGACAUGU	569	21	3
hsa-miR-3934	UCAGGUGUGGAAACUGAGGCAG	570	22	3
hsa-miR-3942-3p	UUUCAGAUAACAGUAUUACAU	414	21	3
hsa-miR-3942-5p	AAGCAAUACUGUUACCUGAAAU	571	22	3
hsa-miR-4523	GACCGAGAGGGCCUCGGCUGU	572	21	3
hsa-miR-4640-5p	UGGGCCAGGGAGCAGCUGGUGGG	573	23	3
hsa-miR-4671-5p	ACCGAAGACUGUGCGCUAAUCU	574	22	3
hsa-miR-4709-5p	ACAACAGUGACUUGCUCUCCAA	575	22	3
hsa-miR-4731-3p	CACACAAGUGGCCCCCAACACU	425	22	3
hsa-miR-4731-5p	UGCUGGGGGCCACAUGAGUGUG	576	22	3
hsa-miR-4762-5p	CCAAAUUCUUGAUCAGAAGCCU	577	21	3
hsa-miR-5010-5p	AGGGGGAUGGCAGAGCAAAAUU	578	22	3
hsa-miR-502-5p	AUCCUUGCUAUCUGGGUGCUA	579	21	3
hsa-miR-548d-5p	AAAAGUAAUUGUGUUUUUGCC	580	22	3

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-548i	AAAAGUAAUUGCGGAUUUUGCC	581	22	3
hsa-miR-548j	AAAAGUAAUUGCGGUCUUUGGU	582	22	3
hsa-miR-5587-3p	GCCCCGGGCAGUGUGAUCAUC	284	21	3
hsa-miR-1225-3p	UGAGCCCCUGUGCCGCCCCCAG	369	22	2
hsa-miR-1227	CGUGCCACCCUUUCCCCAG	583	20	2
hsa-miR-1252	AGAAGGAAAUUGAAUUCAUUUA	371	22	2
hsa-miR-1280	UCCCACCGCUGCCACCC	584	17	2
hsa-miR-1288	UGGACUGCCCUGAUCUGGAGA	585	21	2
hsa-miR-1303	UUUAGAGACGGGGUCUUGCUCU	586	22	2
hsa-miR-1306-3p	ACGUUGGCUCUGGUGGUG	376	18	2
hsa-miR-139-5p	UCUACAGUGCACGUGUCUCCAG	587	22	2
hsa-miR-149-3p	AGGGAGGGACGGGGCUGUGC	588	21	2
hsa-miR-16-1-3p	CCAGUAUUAACUGUGCUGCUGA	589	22	2
hsa-miR-1909-5p	UGAGUGCCGGUGCCUGCCCUG	590	21	2
hsa-miR-224-5p	CAAGUCACUAGUGGUUCCGUU	591	21	2
hsa-miR-2276	UCUGCAAGUGUCAGAGGCGAGG	592	22	2
hsa-miR-2355-3p	AUUGUCCUUGCUGUUUGGAGAU	468	22	2
hsa-miR-2964a-3p	AGAAUUGCGUUUGGACAAUCAGU	392	23	2
hsa-miR-29c-5p	UGACCGAUUUCUCCUGGUGUUC	594	22	2
hsa-miR-3074-3p	GAUAUCAGCUCAGUAGGCACCG	595	22	2
hsa-miR-3120-3p	CACAGCAAGUGUAGACAGGCA	596	21	2
hsa-miR-3130-5p	UACCCAGUCUCCGGUGCAGCC	396	21	2
hsa-miR-3140-3p	AGCUUUUGGGAAUUCAGGUAGU	597	22	2
hsa-miR-3155a	CCAGGCUCUGCAGUGGGAACU	398	21	2
hsa-miR-3163	UAUAAAAUGAGGGCAGUAAGAC	598	22	2
hsa-miR-3167	AGGAUUUCAGAAAUACUGGUGU	599	22	2
hsa-miR-363-5p	CGGGUGGAUCACGAUGCAAUUU	600	22	2
hsa-miR-3676-3p	CCGUGUUUCCCCACGCUUU	408	20	2
hsa-miR-378g	ACUGGGCUUGGAGUCAGAAG	411	20	2
hsa-miR-4467	UGGCGGCGGUAGUUAUGGGCUU	360	22	2
hsa-miR-4498	UGGGCUGGCAGGGCAAGUGCUG	601	22	2
hsa-miR-4654	UGUGGGAUCUGGAGGCAUCUGG	420	22	2
hsa-miR-4659a-3p	UUUCUUCUUAGACAUGGCAACG	603	22	2
hsa-miR-4662a-5p	UUAGCCAAUUGUCCAUCUUUAG	604	22	2
hsa-miR-4683	UGGAGAUCCAGUGCUCGCCCGAU	605	23	2

EP 3 033 418 B1

(continued)

5	CELLS - CTX0E03 07EI				
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
10	hsa-miR-4738-3p	UGAAACUGGAGCGCCUGGAGGA	606	22	2
	hsa-miR-4746-3p	AGCGGUGCUCCUGCGGGCCGA	607	21	2
	hsa-miR-4748	GAGGUUUGGGGAGGAUUUGCU	608	21	2
	hsa-miR-4792	CGGUGAGCGCUCGCUGGC	363	18	2
	hsa-miR-491-5p	AGUGGGGAACCCUCCAUGAGG	429	22	2
15	hsa-miR-5000-3p	UCAGGACACUUCUGAACUUGGA	609	22	2
	hsa-miR-503	UAGCAGCGGGAACAGUUCUGCAG	610	23	2
	hsa-miR-5189	UCUGGGCACAGGCGGAUGGACAGG	611	24	2
	hsa-miR-548aq-3p	CAAAAACUGCAAUUACUUUUGC	612	22	2
	hsa-miR-548av-3p	AAAACUGCAGUUACUUUUGC	613	20	2
20	hsa-miR-5584-5p	CAGGGAAAUGGGAAGAACUAGA	332	22	2
	hsa-miR-5690	UCAGCUACUACCUCUAUUAGG	435	21	2
	hsa-miR-573	CUGAAGUGAUGUGUAACUGAUCAG	305	24	2
	hsa-miR-597	UGUGUCACUCGAUGACCACUGU	614	22	2
	hsa-miR-622	ACAGUCUGCUGAGGUUGGAGC	615	21	2
30	hsa-miR-636	UGUGCUUGCUCGUCCCCGCCGCA	616	23	2
	hsa-miR-1193	GGGAUGGUAGACCGGUGACGUGC	617	23	1
	hsa-miR-1224-3p	CCCCACCUCUCUCUCCUCAG	618	21	1
	hsa-miR-122-5p	UGGAGUGUGACAAUGGUGUUUG	720	22	1
	hsa-miR-1228-5p	GUGGGCGGGGGCAGGUGUGUG	620	21	1
35	hsa-miR-1244	AAGUAGUUGGUUUGUAUGAGAUGGUU	340	26	1
	hsa-miR-1247-5p	ACCCGUCCCGUUCGUCCCCGGA	621	22	1
	hsa-miR-1255b-5p	CGGAUGAGCAAAGAAAGUGGUU	622	22	1
	hsa-miR-1269b	CUGGACUGAGCCAUGCUCUGG	623	22	1
	hsa-miR-1272	GAUGAUGAUGGCAGCAAUUCUGAAA	624	26	1
45	hsa-miR-1273c	GGCGACAAAACGAGACCCUGUC	625	22	1
	hsa-miR-1273e	UUGCUUGAACCCAGGAAGUGGA	342	22	1
	hsa-miR-1282	UCGUUUGCCUUUUUCUGCUU	626	20	1
	hsa-miR-1290	UGGAUUUUUGGAUCAGGGA	375	19	1
	hsa-miR-1294	UGUGAGGUUGGCAUUGUUGUCU	627	22	1
50	hsa-miR-1306-5p	CCACCUCUUUUUGCAAACGUCCA	628	22	1
	hsa-miR-1321	CAGGGAGGUGAAUGUGAU	377	18	1
	hsa-miR-135a-5p	UAUGGCUUUUUUAUCCUAUGUGA	629	23	1
	hsa-miR-137	UUAUUGCUUAAGAAUACGCGUAG	630	23	1
	hsa-miR-142-5p	CAUAAAGUAGAAAGCACUACU	631	21	1

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-143-5p	GGUGCAGUGCUGCAUCUCUGGU	632	22	1
hsa-miR-15a-3p	CAGGCCAUAUUGUGCUGCCUCA	633	22	1
hsa-miR-186-3p	GCCCCAAGGUGAAUUUUUUUGGG	382	22	1
hsa-miR-192-3p	CUGCCAAUCCAUAAGGUCACAG	634	22	1
hsa-miR-19b-1-5p	AGUUUUGCAGGUUUGCAUCCAGC	387	23	1
hsa-miR-200a-3p	UAACACUGUCUGGUAACGAUGU	635	22	1
hsa-miR-204-3p	GCUGGGAAGGCAAAGGGACGU	636	21	1
hsa-miR-214-3p	ACAGCAGGCACAGACAGGCAGU	637	22	1
hsa-miR-29a-5p	ACUGAUUUCUUUUGGUGUUCAG	393	22	1
hsa-miR-3064-5p	UCUGGCUGUUGUGGUGUGCAA	638	21	1
hsa-miR-3116	UGCCUGGAACAUAGUAGGGACU	639	22	1
hsa-miR-3125	UAGAGGAAGCUGUGGAGAGA	640	20	1
hsa-miR-3127-3p	UCCCCUUCUGCAGGCCUGCUGG	641	22	1
hsa-miR-3130-3p	GCUGCACCGGAGACUGGGUAA	395	21	1
hsa-miR-3140-5p	ACCUGAAUUACCAAAAGCUUU	397	21	1
hsa-miR-3157-5p	UUCAGCCAGGCUAGUGCAGUCU	642	22	1
hsa-miR-3179	AGAAGGGGUGAAAUUUAAACGU	643	22	1
hsa-miR-3181	AUCGGGCCCUCGGCGCCGG	644	19	1
hsa-miR-3187-5p	CCUGGGCAGCGUGUGGCUGAAGG	645	23	1
hsa-miR-3190-5p	UCUGGCCAGCUACGUCCCCA	646	20	1
hsa-miR-3198	GUGGAGUCCUGGGGAUUGGAGA	647	22	1
hsa-miR-320b	AAAAGCUGGGUUGAGAGGGCAA	648	22	1
hsa-miR-323b-5p	AGGUUGUCCGUGGUGAGUUCGCA	401	23	1
hsa-miR-3591-5p	UUUAGUGUGAUAAUGGCGUUUGA	649	23	1
hsa-miR-3619-5p	UCAGCAGGCAGGCUGGUGCAGC	650	22	1
hsa-miR-3659	UGAGUGUUGUCUACGAGGGCA	651	21	1
hsa-miR-3674	AUUGUAGAACCUAAGAUUGGCC	652	22	1
hsa-miR-3679-3p	CUUCCCCCAGUAAUCUUCAUC	653	22	1
hsa-miR-375	UUUGUUCGUUCGGCUCGCGUGA	654	22	1
hsa-miR-378b	ACUGGACU UGGAGGCAGAA	655	19	1
hsa-miR-3908	GAGCAAUGUAGGUAGACUGUUU	656	22	1
hsa-miR-3911	UGUGUGGAUCCUGGAGGAGGCA	657	22	1
hsa-miR-3913-5p	UUUGGGACUGAUCUUGAUGUCU	658	22	1
hsa-miR-3917	GCUCGGACUGAGCAGGUGGG	659	20	1
hsa-miR-3944-3p	UUCGGGCUGGCCUGCUGCUCGG	660	23	1

EP 3 033 418 B1

(continued)

	CELLS - CTX0E03 07EI				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-429	UAAUACUGUCUGGUAAAACCGU	661	22	1
	hsa-miR-4421	ACCUGUCUGUGGAAAGGAGCUA	718	22	1
10	hsa-miR-4443	UUGGAGGCGUGGGUUUU	663	17	1
	hsa-miR-4459	CCAGGAGGCGGAGGAGGUGGAG	664	22	1
	hsa-miR-4473	CUAGUGCUCUCCGUUACAAGUA	665	22	1
	hsa-miR-4479	CGCGCGGCCGUGCUCGGAGCAG	666	22	1
15	hsa-miR-4497	CUCCGGGACGGCUGGGC	232	17	1
	hsa-miR-4504	UGUGACAAUAGAGAUGAACAU	667	22	1
	hsa-miR-4520b-3p	UUUGGACAGAAAACACGCAGGU	668	22	1
20	hsa-miR-452-5p	AACUGUUUGCAGAGGAAACUGA	669	22	1
	hsa-miR-4636	AACUCGUGUUCAAAGCCUUUAG	670	22	1
	hsa-miR-4659b-3p	UUUCUUCUAGACAUGGCAGCU	671	22	1
25	hsa-miR-4664-3p	CUUCCGGUCUGUGAGCCCCGUC	672	22	1
	hsa-miR-4665-5p	CUGGGGGACGCGUGAGCGCGAGC	673	23	1
	hsa-miR-4666a-5p	AUACAUGUCAGAUUGUAUGCC	674	21	1
	hsa-miR-4673	UCCAGGCAGGAGCCGGACUGGA	422	22	1
30	hsa-miR-4681	AACGGGAAUGCAGGCUGUAUCU	675	22	1
	hsa-miR-4682	UCUGAGUUCUGGAGCCUGGUCU	676	23	1
	hsa-miR-4690-5p	GAGCAGGCGAGGCUGGGCUGAA	677	22	1
35	hsa-miR-4699-5p	AGAAGAUUGCAGAGUAAGUUC	678	22	1
	hsa-miR-4700-3p	CACAGGACUGACUCCUCACCCCAGUG	424	26	1
	hsa-miR-4706	AGCGGGGAGGAAGUGGGCGCUGCUU	679	25	1
40	hsa-miR-4721	UGAGGGCUCCAGGUGACGGUGG	680	22	1
	hsa-miR-4728-3p	CAUGCUGACCUCUCCUGCCCCAG	681	25	1
	hsa-miR-4742-5p	UCAGGCAAAGGGAUUAUUACAGA	682	23	1
	hsa-miR-4747-3p	AAGGCCCGGGCUUCCUCCAG	683	22	1
45	hsa-miR-4749-5p	UGCGGGGACAGGCCAGGGCAUC	684	22	1
	hsa-miR-4755-3p	AGCCAGGCUCUGAAGGGAAAGU	685	22	1
	hsa-miR-4763-5p	CGCCUGCCCAGCCUCCUGCU	686	21	1
50	hsa-miR-4766-3p	AUAGCAAUUGCUCUUUUGGAA	687	21	1
	hsa-miR-4781-3p	AAUGUUGGAAUCCUCGCUAGAG	688	22	1
	hsa-miR-4793-3p	UCUGCACUGUGAGUUGGCUGGCU	689	23	1
	hsa-miR-488-3p	UUGAAAGGCUAUUUCUUGGUC	690	21	1
55	hsa-miR-4999-5p	UGCUGUAUUGUCAGGUAGUGA	691	21	1
	hsa-miR-5001-5p	AGGGCUGGACUCAGCGCGGAGCU	692	24	1

EP 3 033 418 B1

(continued)

CELLS - CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-5002-5p	AAUUUGGUUUCUGAGGCACUUAGU	693	24	1
hsa-miR-5004-5p	UGAGGACAGGGCAAAUUCACGA	694	22	1
hsa-miR-5006-3p	UUUCCCUUUCUCCUGGCAG	695	21	1
hsa-miR-5088	CAGGGCUCAGGGAUUGGAUGGAG	696	23	1
hsa-miR-544a	AUUCUGCAUUUUUAGCAAGUUC	697	22	1
hsa-miR-548a1	AACGGCAAUGACUUUUGUACCA	698	22	1
hsa-miR-548aq-5p	GAAAGUAAUUGCUGUUUUUGCC	699	22	1
hsa-miR-548at-5p	AAAAGUAAUUGCGUUUUUGGCU	700	22	1
hsa-miR-548au-5p	AAAAGUAAUUGCGUUUUUGC	701	21	1
hsa-miR-548b-3p	CAAGAACCUCAGUUGCUGUUUGU	702	22	1
hsa-miR-556-3p	AUAUUACCAUAGCUCAUCUUU	703	22	1
hsa-miR-5582-3p	UAAAACUUUAAGUGUGCCUAGG	704	22	1
hsa-miR-5586-3p	CAGAGUGACAAGCUGGUUAAAG	705	22	1
hsa-miR-5588-5p	ACUGGCAUUAGUGGGACUUUU	706	21	1
hsa-miR-5683	UACAGAUGCAGAUUCUCUGACUUC	707	24	1
hsa-miR-5696	CUCAUUUAAGUAGUCUGAUGCC	708	22	1
hsa-miR-5701	UUAUUGUCACGUUCUGAUU	709	19	1
hsa-miR-5706	UUCUGGAUACAUGCUGAAGCU	710	22	1
hsa-miR-592	UUGUGUCAUAUGCGAUGAUGU	711	22	1
hsa-miR-603	CACACACUGCAAUUACUUUUGC	712	22	1
hsa-miR-624-3p	CACAAGGUAUUGGUUUACCU	713	21	1
hsa-miR-885-5p	UCCAUUACACUACCCUGCCUCU	714	22	1
hsa-miR-933	UGUGCGCAGGGAGACCUCUCCC	715	22	1

Table 7: Microvesicles EI

MICROVESICLES CTX0E0307EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-1246	AAUGGAUUUUUGGAGCAGG	21	19	32723
hsa-miR-4492	GGGGCUGGGCGCGCGCC	34	17	16225
hsa-miR-4488	AGGGGGCGGGCUCGCGCG	61	18	12878
hsa-miR-4532	CCCCGGGAGCCCGGCG	23	17	6746
hsa-miR-4508	GCGGGGCUGGGCGCGCG	135	17	531
hsa-miR-4516	GGGAGAAGGGUCGGGGC	110	17	500
hsa-miR-3676-5p	AGGAGAUCCUGGGUU	280	15	357

EP 3 033 418 B1

(continued)

	MICROVESICLES CTX0E0307EI				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-4485	UAACGGCCGCGGUACCCUAA	67	20	44
	hsa-miR-4497	CUCCGGGACGGCUGGGC	232	17	43
10	hsa-miR-21-5p	UAGCUUAUCAGACUGAUGUUGA	9	22	33
	hsa-miR-3195	CGCGCCGGGCCCCGGGUU	716	17	28
	hsa-miR-3648	AGCCGCGGGGAUCGCCGAGGG	259	21	26
15	hsa-miR-663b	GGUGGCCCCGGCCGUGCCUGAGG	180	22	24
	hsa-miR-3656	GGCGGGUGCGGGGGUGG	251	17	19
	hsa-miR-3687	CCCGGACAGGCGUUCGUGCGACGU	190	24	19
	hsa-miR-4466	GGGUGCGGGCCGGCGGGG	264	18	19
20	hsa-miR-4792	CGGUGAGCGCUCGCUGGC	363	18	19
	hsa-miR-99b-5p	CACCCGUAGAACCGACCUUGCG	4	22	18
	hsa-let-7a-5p	UGAGGUAGUAGGUUGUAUAGUU	1	22	15
25	hsa-miR-1290	UGGAUUUUUGGAUCAGGGA	375	19	7
	hsa-miR-1291	UGGCCCUGACUGAAGACCAGCAGU	294	24	7
	hsa-miR-182-5p	UUUGGCAAUGGUAGAACUCACACU	16	24	7
	hsa-miR-5096	GUUUCACCAUGUUGGUCAGGC	220	21	7
30	hsa-miR-1273f	GGAGAUGGAGGUUGCAGUG	292	19	5
	hsa-miR-26a-5p	UUCAAGUAAUCCAGGAUAGGCU	12	22	5
	hsa-miR-4284	GGGCUCACAUCACCCCAU	191	18	5
35	hsa-miR-92b-3p	UAUUGCACUCGUCCCCGGCCUCC	13	22	5
	hsa-let-7b-5p	UGAGGUAGUAGGUUGUGUGGUU	28	22	4
	hsa-let-7c	UGAGGUAGUAGGUUGUAUGGUU	17	22	4
	hsa-let-7f-5p	UGAGGUAGUAGAUUGUAUAGUU	11	22	4
40	hsa-miR-100-5p	AACCCGUAGAUCCGAACUUGUG	3	22	4
	hsa-miR-1248	ACCUUCUUGUAUAAGCACUGUGCUGAAA	269	27	4
	hsa-miR-1973	ACCGUGCAAAGGUAGCAUA	171	19	4
45	hsa-miR-21-3p	CAACACCAGUCGAUGGGCUGU	20	21	4
	hsa-miR-3654	GACUGGACAAGCUGAGGAA	325	19	4
	hsa-miR-92a-3p	UAUUGCACUUGUCCCGGCCUGU	7	22	4
50	hsa-miR-1273g-3p	ACCACUGCACUCCAGCCUGAG	210	21	3
	hsa-miR-23b-3p	AUCACAUUGCCAGGGAUUACC	59	21	3
	hsa-miR-3609	CAAAGUGAUGAGUAAUACUGGCUG	216	24	3
	hsa-miR-3615	UCUCUCGGCUCUCCUCGCGGCUC	323	21	3
55	hsa-miR-3653	CUAAGAAGUUGACUGAAG	544	18	3
	hsa-miR-3960	GGCGGCGGCGGAGGCGGGGG	416	20	3

EP 3 033 418 B1

(continued)

	MICROVESICLES CTX0E0307EI				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-4448	GGCUCCUUGGUCUAGGGGUA	231	20	3
	hsa-let-7d-5p	AGAGGUAGUAGGUUGCAUAGUU	92	22	2
10	hsa-miR-16-5p	UAGCAGCACGUAAAUAUUGGCG	29	22	2
	hsa-miR-181a-5p	AACAUUCAACGCUGUCGGUGAGU	15	23	2
	hsa-miR-181b-5p	AACAUUCAUUGCUGUCGGUGGGU	38	23	2
15	hsa-miR-222-3p	AGCUACAUCUGGCUACUGGGU	36	21	2
	hsa-miR-24-3p	UGGCUCAGUUCAGCAGGAACAG	119	22	2
	hsa-miR-3196	CGGGGCGGCAGGGGCCUC	717	18	2
	hsa-miR-4419b	GAGGCUGAAGGAAGAUGG	718	18	2
20	hsa-miR-4461	GAUUGAGACUAGUAGGGCUAGGC	129	23	2
	hsa-miR-4486	GCUGGGCGAGGCUGGCA	719	17	2
	hsa-miR-663a	AGGCGGGGCGCCGCGGGACCGC	365	22	2
25	hsa-miR-9-5p	UCUUUGGUUAUCUAGCUGUAUGA	58	23	2
	hsa-let-7i-3p	CUGCGCAAGCUACUGCCUUGCU	483	22	1
	hsa-let-7i-5p	UGAGGUAGUAGUUUGUGCUGUU	22	22	1
	hsa-miR-1225-5p	GUGGGUACGGCCCAGUGGGGGG	720	22	1
30	hsa-miR-1244	AAGUAGUUGGUUUUGUAUGAGAUGGUU	340	26	1
	hsa-miR-125b-5p	UCCCUGAGACCCUAACUUGUGA	42	22	1
	hsa-miR-1275	GUGGGGGAGAGGCUGUC	162	17	1
35	hsa-miR-1280	UCCCACCGCUGCCACCC	584	17	1
	hsa-miR-134	UGUGACUGGUUGACCAGAGGGG	94	22	1
	hsa-miR-149-5p	UCUGGCUCCGUGUCUUCACUCCC	121	23	1
	hsa-miR-191-5p	CAACGGAAUCCCAAAAGCAGCUG	8	23	1
40	hsa-miR-221-3p	AGCUACAUUGUCUGCUGGGUUUC	79	23	1
	hsa-miR-22-3p	AAGCUGCCAGUUGAAGAACUGU	33	22	1
	hsa-miR-26b-5p	UUCAAGUAAUUCAGGAUAGGU	90	21	1
45	hsa-miR-30c-5p	UGUAAACAUCCUACACUCUCAGC	66	23	1
	hsa-miR-30d-5p	UGUAAACAUCCCCGACUGGAAG	31	22	1
	hsa-miR-3182	GCUUCUGUAGUGUAGUC	721	17	1
50	hsa-miR-320a	AAAAGCUGGGUUGAGAGGGCGA	97	22	1
	hsa-miR-34a-5p	UGGCAGUGUCUUAGCUGGUUGU	101	22	1
	hsa-miR-3607-3p	ACUGUAAACGCUUUCUGAUG	543	20	1
	hsa-miR-361-5p	UUAUCAGAAUCUCCAGGGGUAC	70	22	1
55	hsa-miR-3652	CGGCUGGAGGUGUGAGGA	722	18	1
	hsa-miR-409-3p	GAAUGUUGCUCGGUGAACCCCU	47	22	1

EP 3 033 418 B1

(continued)

MICROVESICLES CTX0E0307EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-423-3p	AGCUCGGUCUGAGGCCCCUCAGU	57	23	1
hsa-miR-423-5p	UGAGGGGCAGAGAGCGAGACUUU	41	23	1
hsa-miR-432-5p	UCUUGGAGUAGGUCAUUGGGUGG	95	23	1
hsa-miR-4417	GGUGGGCUUCCCGGAGGG	175	18	1
hsa-miR-4426	GAAGAUGGACGUACUUU	359	17	1
hsa-miR-4449	CGUCCCGGGGUCGCGGAGGCA	155	22	1
hsa-miR-4800-3p	CAUCCGUCCGUCUGUCCAC	549	19	1
hsa-miR-484	UCAGGCUCAGUCCCCUCCGAU	118	22	1
hsa-miR-486-5p	UCCUGUACUGAGCUGCCCCGAG	5	22	1
hsa-miR-493-3p	UGAAGGUCUACUGUGUGCCAGG	83	22	1
hsa-miR-5095	UUACAGGCGUGAACCACCGCG	723	21	1
hsa-miR-556-3p	AUAUUACCAUUAGCUCAUCUUU	703	22	1
hsa-miR-644b-5p	UGGGCUAAGGGAGAUGAUUGGGUA	724	24	1
hsa-miR-664-5p	ACUGGCUAGGGAAAAUGAUUGGAU	443	24	1
hsa-miR-760	CGGCUCUGGGUCUGUGGGGA	289	20	1
hsa-miR-941	CACCCGGCUGUGUGCACAUGUGC	60	23	1
hsa-miR-98	UGAGGUAGUAAGUUGUAUUGUU	10	22	1
hsa-miR-99a-5p	AACCCGUAGAUCCGAUCUUGUG	52	22	1

Table 8: Exosomes EI

EXOSOMES CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-1246	AAUGGAUUUUUGGAGCAGG	21	19	83958
hsa-miR-4492	GGGGCUGGGCGCGCGCC	34	17	22482
hsa-miR-4488	AGGGGGCGGGCUCGGCG	61	18	20618
hsa-miR-4532	CCCCGGGGAGCCCGGCG	23	17	6419
hsa-miR-4516	GGGAGAAGGGUCGGGGC	110	17	904
hsa-miR-4508	GCGGGGUCUGGGCGCGCG	135	17	723
hsa-miR-3676-5p	AGGAGAUCCUGGGUU	280	15	174
hsa-miR-4485	UAACGGCCGCGGUACCCUAA	67	20	43
hsa-miR-21-5p	UAGCUUAUCAGACUGAUUGUA	9	22	41
hsa-miR-4497	CUCCGGGACGGCUGGGC	232	17	28
hsa-miR-663b	GGUGCCCCGGCCGUGCCUGAGG	180	22	26
hsa-miR-4792	CGGUGAGCGCUCGCUGGC	363	18	24

EP 3 033 418 B1

(continued)

EXOSOMES CTX0E03 07EI				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-4454	GGAUCCGAGUCACGGCACCA	299	20	22
hsa-miR-1291	UGGCCUGACUGAAGACCAGCAGU	294	24	17
hsa-miR-26a-5p	UUCAAGUAAUCCAGGAUAGGCU	12	22	17
hsa-miR-3195	CGCGCCGGGCCCCGGGUU	716	17	17
hsa-let-7a-5p	UGAGGUAGUAGGUUGUAUAGUU	1	22	15
hsa-miR-182-5p	UUUGGCAAUGGUAGAACUCACACU	16	24	15
hsa-miR-99b-5p	CACCCGUAGAACCGACCUUGCG	4	22	15
hsa-miR-5096	GUUUCACCAUGUUGGUCAGGC	220	21	14
hsa-miR-3648	AGCCGCGGGGAUCGCCGAGGG	259	21	13
hsa-miR-3654	GACUGGACAAGCUGAGGAA	325	19	13
hsa-miR-4466	GGGUGCGGGCCGCGGGG	264	18	12
hsa-miR-3687	CCCGGACAGGCGUUCGUGCGACGU	190	24	11
hsa-miR-4284	GGGCUCACAUACCCCAU	191	18	11
hsa-miR-3656	GGCGGGUGCGGGGGUGG	251	17	10
hsa-miR-3609	CAAAGUGAUGAGUAAUACUGGCUG	216	24	8
hsa-miR-644b-5p	UGGGCUAAGGGAGAUGAUUGGGUA	724	24	8
hsa-miR-664-5p	ACUGGCUAGGGAAAAUGAUUGGAU	443	24	8
hsa-miR-92a-3p	UAUUGCACUUGUCCCGGCCUGU	7	22	7
hsa-miR-92b-3p	UAUUGCACUCGUCCCGGCCUCC	13	22	7
hsa-let-7b-5p	UGAGGUAGUAGGUUGUGUGGUU	28	22	6
hsa-let-7f-5p	UGAGGUAGUAGAUUGUAUAGUU	11	22	6
hsa-miR-127-3p	UCGGAUCCGUCUGAGCUUGGCU	14	22	6
hsa-miR-1290	UGGAUUUUUGGAUCAGGGA	375	19	6
hsa-miR-4449	CGUCCCGGGGCU GCGCGAGGCA	155	22	6
hsa-miR-4461	GAUUGAGACUAGUAGGGCUAGGC	129	23	6
hsa-miR-100-5p	AACCCGUAGAUCCGAACUUGUG	3	22	5
hsa-miR-1248	ACCUUCUUGUAUAAGCACUGUGC AAA	269	27	5
hsa-miR-1973	ACCGUGCAAAGGUAGCAUA	171	19	5
hsa-miR-3653	CUAAGAAGUUGACUGAAG	544	18	5
hsa-miR-4417	GGUGGGCUUCCCGGAGGG	175	18	5
hsa-miR-125b-5p	UCCCUGAGACCCUAACUUGUGA	42	22	4
hsa-miR-151a-3p	CUAGACUGAAGCUCCUUGAGG	25	21	4
hsa-miR-16-5p	UAGCAGCACGUAAUUAUUGGCG	29	22	4
hsa-miR-21-3p	CAACACCAGUCGAUGGGCUGU	20	21	4
hsa-miR-23a-3p	AUCACAUUGCCAGGGAUUUCC	55	21	4

EP 3 033 418 B1

(continued)

	EXOSOMES CTX0E03 07EI				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-4419b	GAGGCUGAAGGAAGAUGG	718	18	4
10	hsa-miR-1273f	GGAGAUGGAGGUUGCAGUG	292	19	3
	hsa-miR-1273g-3p	ACCACUGCACUCCAGCCUGAG	210	21	3
	hsa-miR-181b-5p	AACAUUCAUUGCUGUCGGUGGGU	38	23	3
	hsa-miR-221-3p	AGCUACAUUGCUGCUGGGUUUC	79	23	3
15	hsa-miR-3615	UCUCUCGGCUCUCGCGGCUC	323	21	3
	hsa-miR-9-5p	UCUUUGGUUAUCUAGCUGUAUGA	58	23	3
	hsa-let-7c	UGAGGUAGUAGGUUGUAUGGUU	17	22	2
20	hsa-let-7e-5p	UGAGGUAGGAGGUUGUAUAGUU	27	22	2
	hsa-let-7i-5p	UGAGGUAGUAGUUUGUCUGUU	22	22	2
	hsa-miR-103a-3p	AGCAGCAUUGUACAGGGCUAUGA	62	23	2
	hsa-miR-106b-5p	UAAAGUGCUGACAGUGCAGAU	170	21	2
25	hsa-miR-1273e	UUGCUUGAACCCAGGAAGUGGA	342	22	2
	hsa-miR-221-5p	ACCUGGCAUACAAUGUAGAUUU	39	22	2
	hsa-miR-222-3p	AGCUACAUCUGGCUACUGGGU	36	21	2
30	hsa-miR-30d-5p	UGUAAACAUCCCCGACUGGAAG	31	22	2
	hsa-miR-3960	GGCGGCGGCGGAGGCGGGGG	416	20	2
	hsa-let-7d-3p	CUAUACGACCUGCUGCCUUUCU	92	22	1
	hsa-let-7d-5p	AGAGGUAGUAGGUUGCAUAGUU	53	22	1
35	hsa-let-7g-5p	UGAGGUAGUAGUUUGUACAGUU	43	22	1
	hsa-let-7i-3p	CUGCGCAAGCUACUGCCUUGCU	483	22	1
	hsa-miR-10a-5p	UACCCUGUAGAUCCGAAUUUGUG	2	23	1
40	hsa-miR-1181	CCGUCGCCGCCACCCGAGCCG	725	21	1
	hsa-miR-1225-3p	UGAGCCCCUGUGCCGCCCCAG	369	22	1
	hsa-miR-1244	AAGUAGUUGGUUUGUAUGAGAUGGUU	340	26	1
	hsa-miR-125a-5p	UCCCUGAGACCCUUUAACCGUGA	35	24	1
45	hsa-miR-1296	UUAGGGCCCUGGCUCCAUCUCC	271	22	1
	hsa-miR-1307-5p	UCGACCGGACCUCGACCGGCU	91	21	1
	hsa-miR-146b-5p	UGAGAACUGAAUCCAUAAGGCU	19	22	1
50	hsa-miR-149-5p	UCUGGCUCCGUGUCUUCACUCCC	121	23	1
	hsa-miR-151a-5p	UCGAGGAGCUCACAGUCUAGU	37	21	1
	hsa-miR-15b-5p	UAGCAGCACAUCAUGGUUUACA	78	22	1
	hsa-miR-181a-2-3p	ACCACUGACCGUUGACUGUACC	102	22	1
55	hsa-miR-181a-5p	AACAUUCAACGCUGUCGGUGAGU	15	23	1
	hsa-miR-191-5p	CAACGGAAUCCCAAAGCAGCUG	8	23	1

EP 3 033 418 B1

(continued)

	EXOSOMES CTX0E03 07EI				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-198	GGUCCAGAGGGGAGAUAGGUUC	726	22	1
10	hsa-miR-204-5p	UCCCCUUUGUCAUCCU AUGCCU	89	22	1
	hsa-miR-20a-5p	UAAAGUGCUUUAUAGUGCAGGUAG	146	23	1
	hsa-miR-219-5p	UGAUUGUCCAAACGCAAUUCU	527	21	1
	hsa-miR-22-3p	AAGCUGCCAGUUGAAGAACUGU	33	22	1
15	hsa-miR-23b-3p	AUCACAUUGCCAGGGAUUACC	59	21	1
	hsa-miR-26b-3p	CCUGUUCUCCAUUACUUGGCUC	391	22	1
	hsa-miR-299-5p	UGGUUUACCGUCCACAUACA	319	22	1
20	hsa-miR-29a-3p	UAGCACCAUCUGAAAUCGGUUA	106	22	1
	hsa-miR-30e-3p	CUUUCAGUCGGAUGUUUACAGC	71	22	1
	hsa-miR-31-3p	UGCUAUGCCAACAUAUUGCCAU	172	22	1
	hsa-miR-3198	GUGGAGUCCUGGGGAUUGGAGA	647	22	1
25	hsa-miR-323a-3p	CACAUUACACGGUCGACCUCU	158	21	1
	hsa-miR-342-3p	UCUCACACAGAAAUCGCACCCGU	81	23	1
	hsa-miR-3607-3p	ACUGUAAACGCUUUCUGAUG	543	20	1
30	hsa-miR-3651	CAUAGCCCCGGUCGCGUGUACAUGA	727	24	1
	hsa-miR-378a-3p	ACUGGACUUGGAGUCAGAAGG	65	21	1
	hsa-miR-379-5p	UGGUAGACUAUGGAACGUAGG	18	21	1
	hsa-miR-423-3p	AGCUCGGUCUGAGGCCCCUCAGU	57	23	1
35	hsa-miR-423-5p	UGAGGGGCAGAGAGCGAGACUUU	41	23	1
	hsa-miR-425-5p	AAUGACACGAUCACUCCCGUUGA	111	23	1
	hsa-miR-4258	CCCCGCCACCGCCUUGG	728	17	1
40	hsa-miR-4426	GAAGAUGGACGUACUUU	359	17	1
	hsa-miR-4443	UUGGAGGCGUGGGUUUU	663	17	1
	hsa-miR-4448	GGCUCCUUGGUCUAGGGGUA	231	20	1
	hsa-miR-4697-3p	UGUCAGUGACUCCUGCCCCUUGGU	729	24	1
45	hsa-miR-4700-3p	CACAGGACUGACUCCUCACCCAGUG	424	26	1
	hsa-miR-4700-5p	UCUGGGGAUGAGGACAGUGUGU	730	22	1
	hsa-miR-4797-3p	UCUCAGUAAGUGGCACUCUGU	731	21	1
50	hsa-miR-484	UCAGGCUCAGUCCCCUCCGAU	118	22	1
	hsa-miR-486-5p	UCCUGUACUGAGCUGCCCCGAG	5	22	1
	hsa-miR-494	UGAAACAUACACGGGAAACCUC	240	22	1
	hsa-miR-500a-5p	UAAUCCUUGCUACCUGGGUGAGA	303	23	1
55	hsa-miR-644b-3p	UUCAUUUGCCUCCAGCCUACA	442	22	1
	hsa-miR-663a	AGGCGGGGCGCCGCGGGACCGC	365	22	1

Table 9: Microvesicles EH

5	MICROVESICLES CTX0E03 07EH				
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
10	hsa-miR-1246	AAUGGAUUUUUGGAGCAGG	21	19	78791
	hsa-miR-4492	GGGGCUGGGCGCGCGCC	34	17	6012
	hsa-miR-4532	CCCCGGGGAGCCCGGCG	23	17	3410
	hsa-miR-4488	AGGGGGCGGGCUCCGGCG	61	18	1737
15	hsa-miR-4485	UAACGGCCGCGGUACCCUAA	67	20	319
	hsa-miR-4508	GCGGGGCUGGGCGCGCG	135	17	221
	hsa-miR-4516	GGGAGAAGGGUCGGGGC	110	17	114
	hsa-miR-4466	GGGUGCGGGCCGGCGGGG	264	18	61
20	hsa-miR-4497	CUCCGGGACGGCUGGGC	232	17	51
	hsa-miR-3195	CGCGCCGGGCCCCGGGUU	716	17	41
	hsa-miR-1973	ACCGUGCAAAGGUAGCAUA	171	19	30
	hsa-miR-21-5p	UAGCUUAUCAGACUGAUGUUGA	9	22	22
25	hsa-miR-4284	GGGCUCACAUCACCCCAU	191	18	20
	hsa-miR-4792	CGGUGAGCGCUCGCGUGGC	363	18	12
	hsa-miR-92a-3p	UAUUGCACUUGUCCCGGCCUGU	7	22	12
	hsa-miR-1291	UGGCCUGACUGAAGACCAGCAGU	294	24	11
30	hsa-miR-3676-5p	AGGAGAUCCUGGGUU	280	15	10
	hsa-miR-100-5p	AACCCGUAGAUCGGAACUUGUG	3	22	8
	hsa-miR-3656	GGCGGGUGCGGGGGUGG	251	17	8
	hsa-miR-663b	GGUGGCCCGGCCGUGCCUGAGG	180	22	8
35	hsa-let-7a-5p	UGAGGUAGUAGGUUGUAUAGUU	1	22	7
	hsa-miR-1290	UGGAUUUUUGGAUCAGGGA	375	19	7
	hsa-miR-3687	CCCGGACAGGCGUUCGUGCGACGU	190	24	7
	hsa-miR-4461	GAUUGAGACUAGUAGGGCUAGGC	52	23	6
40	hsa-miR-664-5p	ACUGGCUAGGGAAAUGAUUGGAU	91	24	6
	hsa-miR-92b-3p	UAUUGCACUCGUCCCGGCCUCC	13	22	6
	hsa-miR-125b-5p	UCCCUGAGACCCUAACUUGUGA	42	22	5
	hsa-miR-3653	CUAAGAAGUUGACUGAAG	544	18	5
45	hsa-let-7f-5p	UGAGGUAGUAGAUUGUAUAGUU	11	22	4
	hsa-miR-16-5p	UAGCAGCACGUAAAUAUUGGCG	29	22	4
	hsa-miR-181a-5p	AACAUUCAACGCUGUCGGUGAGU	15	23	4
	hsa-miR-3609	CAAAGUGAUGAGUAAUACUGGCUG	216	24	4
50	hsa-miR-9-5p	UCUUUGGUUAUCUAGCUGUAUGA	58	23	4
	hsa-let-7c	UGAGGUAGUAGGUUGUAUGGUU	17	22	3
	hsa-miR-1244	AAGUAGUUGGUUUUGUAUGAGAUGGUU	59	26	3

EP 3 033 418 B1

(continued)

	MICROVESICLES CTX0E03 07EH				
5	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
	hsa-miR-127-3p	UCGGAUCCGUCUGAGCUUGGCU	14	22	3
	hsa-miR-181b-5p	AACAUUCAUUGCUGUCGGUGGGU	38	23	3
10	hsa-miR-21-3p	CAACACCAGUCGAUGGGCUGU	20	21	3
	hsa-miR-26a-5p	UUCAAGUAAUCCAGGAUAGGCU	12	22	3
	hsa-miR-30c-5p	UGUAAACAUCCUACACUCUCAGC	66	23	3
15	hsa-miR-3960	GGCGGCGGCGGAGGCGGGGG	416	20	3
	hsa-miR-485-3p	GUCAUACACGGCUCUCCUCUCU	153	22	3
	hsa-let-7b-5p	UGAGGUAGUAGGUUGUGUGGUU	28	22	2
	hsa-let-7g-5p	UGAGGUAGUAGUUUGUACAGUU	43	22	2
20	hsa-miR-1273f	GGAGAUGGAGGUUGCAGUG	292	19	2
	hsa-miR-151a-3p	CUAGACUGAAGCUCCUUGAGG	25	21	2
	hsa-miR-182-5p	UUUGGCAAUGGUAGAACUCACACU	16	24	2
25	hsa-miR-191-5p	CAACGGAAUCCCAAAGCAGCUG	8	23	2
	hsa-miR-197-3p	UUCACCACCUUCUCCACCCAGC	122	22	2
	hsa-miR-423-5p	UGAGGGGCAGAGAGCGAGACUUU	41	23	2
	hsa-miR-4468	AGAGCAGAAGGAUGAGAU	732	18	2
30	hsa-miR-644b-5p	UGGGCUAAGGGAGAUGAUUGGGUA	724	24	2
	hsa-miR-93-5p	CAAAGUGCUGUUCGUGCAGGUAG	116	23	2
	hsa-let-7d-5p	AGAGGUAGUAGGUUGCAUAGUU	92	22	1
35	hsa-miR-1225-3p	UGAGCCCCUGUGCCGCCCCAG	369	22	1
	hsa-miR-1254	AGCCUGGAAGCUGGAGCCUGCAGU	270	24	1
	hsa-miR-1273g-3p	ACCACUGCACUCCAGCCUGAG	210	21	1
40	hsa-miR-1275	GUGGGGGAGAGGCUGUC	162	17	1
	hsa-miR-1296	UUAGGGCCCUGGCUCCAUCUCC	271	22	1
	hsa-miR-1307-5p	UCGACCGGACCUCGACCGGCU	91	21	1
	hsa-miR-134	UGUGACUGGUUGACCAGAGGGG	94	22	1
45	hsa-miR-15b-5p	UAGCAGCACAUCAUGGUUUACA	78	22	1
	hsa-miR-17-5p	CAAAGUGCUIACAGUGCAGGUAG	145	23	1
	hsa-miR-1972	UCAGGCCAGGCACAGUGGCUCA	733	22	1
50	hsa-miR-22-3p	AAGCUGCCAGUUGAAGAACUGU	33	22	1
	hsa-miR-25-3p	CAUUGCACUUGUCUCGGUCUGA	63	22	1
	hsa-miR-27b-3p	UUCACAGUGGCUAAGUUCUGC	6	21	1
	hsa-miR-3065-5p	UCAACAAAAUCACUGAUGCUGGA	226	23	1
55	hsa-miR-30d-5p	UGUAAACAUCCCCGACUGGAAG	31	22	1
	hsa-miR-320a	AAAAGCUGGGUUGAGAGGGCGA	97	22	1

EP 3 033 418 B1

(continued)

MICROVESICLES CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-342-3p	UCUCACACAGAAAUCGCACCCGU	81	23	1
hsa-miR-3648	AGCCGCGGGGAUCGCCGAGGG	259	21	1
hsa-miR-3652	CGGCUGGAGGUGUGAGGA	722	18	1
hsa-miR-376c	AACAUAGAGGAAAUUCCACGU	185	21	1
hsa-miR-378a-3p	ACUGGACUUGGAGUCAGAAGG	65	21	1
hsa-miR-409-3p	GAAUGUUGCUCGGUGAACCCCU	47	22	1
hsa-miR-433	AUCAUGAUGGGCUCCUCGGUGU	174	22	1
hsa-miR-4417	GGUGGGCUUCCCGGAGGG	175	18	1
hsa-miR-4448	GGCUCCUUGGUCUAGGGGUA	231	20	1
hsa-miR-4454	GGAUCCGAGUCACGGCACCA	299	20	1
hsa-miR-454-3p	UAGUGCAAUAUUGCUUAUAGGGU	169	23	1
hsa-miR-4800-3p	CAUCCGUCCGUCUGUCCAC	549	19	1
hsa-miR-493-3p	UGAAGGUCUACUGUGGCCAGG	83	22	1
hsa-miR-5095	UUACAGGCGUGAACCAACCGCG	723	21	1
hsa-miR-574-3p	CACGCUCAUGCACACACCCACA	253	22	1
hsa-miR-665	ACCAGGAGGCUGAGGCCCCU	309	20	1
hsa-miR-720	UCUCGCUGGGGCCUCCA	84	17	1
hsa-miR-99a-5p	AACCCGUAGAUCCGAUCUUGUG	52	22	1
hsa-miR-99b-5p	CACCCGUAGAACCGACCUUGCG	4	22	1

Table 10: Exosomes EH

EXOSOMES CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-1246	AAUGGAUUUUUGGAGCAGG	21	19	111092
hsa-miR-4492	GGGGCUGGGCGCGCGCC	34	17	5188
hsa-miR-4532	CCCCGGGGAGCCCGGCG	23	17	3368
hsa-miR-4488	AGGGGGCGGGCUCCGGCG	61	18	1389
hsa-miR-4485	UAACGGCCGCGGUACCCUAA	67	20	386
hsa-miR-4508	GCGGGGCUGGGCGCGCG	135	17	188
hsa-miR-4516	GGGAGAAGGGUCGGGGC	110	17	135
hsa-miR-4497	CUCCGGGACGGCUGGGC	232	17	73
hsa-miR-1973	ACCGUGCAAAGGUAGCAUA	171	19	50
hsa-miR-3195	CGCGCCGGGCCCGGGUU	716	17	48
hsa-miR-4466	GGGUGCGGGCCGGCGGGG	264	18	43

EP 3 033 418 B1

(continued)

5	EXOSOMES CTX0E03 07EH				
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
10	hsa-let-7a-5p	UGAGGUAGUAGGUUGUAUAGUU	1	22	20
	hsa-miR-99b-5p	CACCCGUAGAACCGACCUUGCG	4	22	19
	hsa-miR-21-5p	UAGCUUAUCAGACUGAUGUUGA	9	22	18
	hsa-miR-92a-3p	UAUUGCACUUGUCCCGGCCUGU	7	22	18
	hsa-miR-3676-5p	AGGAGAUCCUGGGUU	280	15	17
15	hsa-miR-4792	CGGUGAGCGCUCGCUGGC	363	18	15
	hsa-miR-664-5p	ACUGGCUAGGGAAAAUGAUUGGAU	443	24	13
	hsa-miR-100-5p	AACCCGUAGAUCCGAACUUGUG	3	22	11
	hsa-miR-1291	UGGCCUGACUGAAGACCAGCAGU	294	24	10
	hsa-miR-16-5p	UAGCAGCACGUAAAUAUUGGCG	29	22	10
20	hsa-miR-4284	GGGCUCACAUCACCCCAU	191	18	10
	hsa-miR-663b	GGUGGCCCGGCCGUGCCUGAGG	180	22	9
	hsa-miR-25-3p	CAUUGCACUUGUCUCGGUCUGA	63	22	8
	hsa-miR-3656	GGCGGGUGCGGGGGUGG	251	17	8
	hsa-miR-181a-5p	AACAUUCAACGCUGUCGGUGAGU	15	23	7
30	hsa-miR-26a-5p	UUCAAGUAAUCCAGGAUAGGCU	12	22	6
	hsa-miR-3654	GACUGGACAAGCUGAGGAA	325	19	6
	hsa-miR-644b-5p	UGGGCUAAGGGAGAUGAUUGGGUA	724	24	6
	hsa-let-7b-5p	UGAGGUAGUAGGUUGUGUGGUU	28	22	5
	hsa-let-7f-5p	UGAGGUAGUAGAUUGUAUAGUU	11	22	5
35	hsa-miR-1290	UGGAUUUUUGGAUCAGGGA	375	19	5
	hsa-miR-4426	GAAGAUGGACGUACUUU	359	17	5
	hsa-miR-5096	GUUUCACCAUGUUGGUCAGGC	220	21	5
	hsa-miR-125b-5p	UCCCUGAGACCCUAACUUGUGA	42	22	4
	hsa-miR-1273f	GGAGAUGGAGGUUGCAGUG	292	19	4
45	hsa-miR-191-5p	CAACGGAAUCCCAAAGCAGCUG	8	23	4
	hsa-miR-22-3p	AAGCUGCCAGUUGAAGAACUGU	33	22	4
	hsa-miR-3609	CAAAGUGAUGAGUAAUACUGGCUG	216	24	4
	hsa-miR-3687	CCCGGACAGGCGUUCGUGCGACGU	190	24	4
	hsa-miR-93-5p	CAAAGUGCUGUUCGUGCAGGUAG	116	23	4
50	hsa-miR-1248	ACCUUCUUGUAUAAGCACUGUGCUAAA	269	27	3
	hsa-miR-1273g-3p	ACCACUGCACUCCAGCCUGAG	210	21	3
	hsa-miR-151a-3p	CUAGACUGAAGCUCCUUGAGG	25	21	3
	hsa-miR-182-5p	UUUGGCAAUGGUAGAACUCACACU	16	24	3
	hsa-miR-221-3p	AGCUACAUUGUCUGCUGGGUUUC	79	23	3

EP 3 033 418 B1

(continued)

5	EXOSOMES CTX0E03 07EH				
	MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
10	hsa-miR-222-3p	AGCUACAUCUGGCUACUGGGU	36	21	3
	hsa-miR-29a-3p	UAGCACCAUCUGAAAUCGGUUA	106	22	3
	hsa-miR-4461	GAUUGAGACUAGUAGGGCUAGGC	129	23	3
	hsa-miR-486-5p	UCCUGUACUGAGCUGCCCCGAG	5	22	3
	hsa-miR-92b-3p	UAUUGCACUCGUCCCGGCCUCC	13	22	3
15	hsa-miR-9-5p	UCUUUGGUUAUCUAGCUGUAUGA	58	23	3
	hsa-miR-98	UGAGGUAGUAAGUUGUAUUGUU	10	22	3
	hsa-let-7d-5p	AGAGGUAGUAGGUUGCAUAGUU	53	22	2
	hsa-miR-134	UGUGACUGGUUGACCAGAGGGG	94	22	2
	hsa-miR-151a-5p	UCGAGGAGCUCACAGUCUAGU	37	21	2
20	hsa-miR-15b-5p	UAGCAGCACAUCAUGGUUUACA	78	22	2
	hsa-miR-30a-5p	UGUAAACAUCCUCGACUGGAAG	30	22	2
	hsa-miR-3124-3p	ACUUUCCUCACUCCCGUGAAGU	734	22	2
	hsa-miR-3653	CUAAGAAGUUGACUGAAG	544	18	2
	hsa-let-7c	UGAGGUAGUAGGUUGUAUGGUU	17	22	1
30	hsa-let-7d-3p	CUAUACGACCUGCUGCCUUUCU	92	22	1
	hsa-let-7g-5p	UGAGGUAGUAGUUUGUACAGUU	43	22	1
	hsa-let-7i-5p	UGAGGUAGUAGUUUGUGCUGUU	22	22	1
	hsa-miR-103a-3p	AGCAGCAUUGUACAGGGCUAUGA	62	23	1
	hsa-miR-106b-5p	UAAAGUGCUGACAGUGCAGAU	170	21	1
35	hsa-miR-1244	AAGUAGUUGGUUUUGUAUGAGAUGGUU	340	26	1
	hsa-miR-128	UCACAGUGAACCGGUCUCUUU	109	21	1
	hsa-miR-1285-3p	UCUGGGCAACAAAGUGAGACCU	464	22	1
	hsa-miR-1307-3p	ACUCGGCGUGGCGUCGGUCGUG	124	22	1
	hsa-miR-140-3p	UACCACAGGGUAGAACCACGG	138	21	1
45	hsa-miR-148b-3p	UCAGUGCAUCACAGAACUUUGU	48	22	1
	hsa-miR-181b-5p	AACAUUCAUUGCUGUCGGUGGGU	38	23	1
	hsa-miR-193a-3p	AACUGGCCUACAAAGUCCCAGU	386	22	1
	hsa-miR-1972	UCAGGCCAGGCACAGUGGCUCA	733	22	1
	hsa-miR-21-3p	CAACACCAGUCGAUGGGCUGU	20	21	1
50	hsa-miR-2277-3p	UGACAGCGCCCUGCCUGGCUC	735	21	1
	hsa-miR-23a-3p	AUCACAUUGCCAGGGAUUUCC	55	21	1
	hsa-miR-23b-3p	AUCACAUUGCCAGGGAUUACC	59	21	1
	hsa-miR-24-3p	UGGCUCAGUUCAGCAGGAACAG	119	22	1
	hsa-miR-27a-3p	UUCACAGUGGCUAAGUUCGCG	46	21	1

EP 3 033 418 B1

(continued)

EXOSOMES CTX0E03 07EH				
MIRNA	MIRNA.SEQUENCE	SEQ ID NO:	MIRNA LENGTH	READ COUNTS
hsa-miR-27b-3p	UUCACAGUGGCUAAGUUCUGC	6	21	1
hsa-miR-299-3p	UAUGUGGGAUGGUAAACCGCUU	182	22	1
hsa-miR-30b-5p	UGUAAACAUCUACACUCAGCU	96	22	1
hsa-miR-30c-5p	UGUAAACAUCUACACUCUCAGC	66	23	1
hsa-miR-31-3p	UGCUAUGCCAACAUAUUGCCAU	172	22	1
hsa-miR-3196	CGGGGCGGCAGGGGCCUC	717	18	1
hsa-miR-3198	GUGGAGUCCUGGGGAUUGGAGA	647	22	1
hsa-miR-320a	AAAAGCUGGGUUGAGAGGGCGA	97	22	1
hsa-miR-329	AACACACCUGGUUAACCUCUUU	214	22	1
hsa-miR-339-5p	UCCUGUCCUCCAGGAGCUCACG	402	23	1
hsa-miR-34a-5p	UGGCAGUGUCUUAGCUGGUUGU	101	22	1
hsa-miR-3607-5p	GCAUGUGAUGAAGCAAAUCAGU	249	22	1
hsa-miR-3648	AGCCGCGGGGAUCGCCGAGGG	259	21	1
hsa-miR-376c	AACAUAGAGGAAAUUCCACGU	185	21	1
hsa-miR-3960	GGCGGCGGCGGAGGCGGGGG	416	20	1
hsa-miR-411-3p	UAUGUAAACACGGUCCACUAACC	482	22	1
hsa-miR-423-3p	AGCUCGGUCUGAGGCCCCUCAGU	57	23	1
hsa-miR-423-5p	UGAGGGGCAGAGAGCGAGACUUU	41	23	1
hsa-miR-4417	GGUGGGCUUCCCGGAGGG	175	18	1
hsa-miR-4444	CUCGAGUUGGAAGAGGCG	418	18	1
hsa-miR-4499	AAGACUGAGAGGAGGGA	736	17	1
hsa-miR-4521	GCUAAGGAAGUCCUGUGCUCAG	233	22	1
hsa-miR-4680-5p	AGAACUCUUGCAGUCUAGAUGU	737	23	1
hsa-miR-4709-5p	ACAACAGUGACUUGCUCUCCAA	575	22	1
hsa-miR-501-3p	AAUGCACCCGGGCAAGGAUUCU	26	22	1
hsa-miR-644b-3p	UUCAUUUGCCUCCAGCCUACA	442	22	1
hsa-miR-654-3p	UAUGUCUGCUGACCAUCACCUU	336	22	1
hsa-miR-9-3p	AUAAAGCUAGAUAAACCGAAAGU	183	22	1
hsa-miR-940	AAGGCAGGGCCCCCGCUCCCC	366	21	1
hsa-miR-99a-5p	AACCCGUAGAUCCGAUCUUGUG	52	22	1

Identification of top ranking coding and non-coding RNAs by GENCODE analysis performed in exosomes, MV and producer cells

[0296]

Table 11: Total number of sequence reads identified by using GENCODE in each tested samples

CTX0E0307EH cells	CTX0E0307EH EXO	CTX0E0307EH MV	CTX0E0307EI cells	CTX0E0307EIE XO	CTX0E0307EI MV
18741941	12678688	10876797	22116110	16311289	835970

[0297] Using GENCODE database analysis of the sequence results, seven putative novel miRNA sequences were identified in exosomes (EXO), microvesicles (MV) and producer cells, as shown in Table 12. (nb CTX0E03 07EI MV reads are misrepresented due to the lower amount of starting material - see Table 11). These data are shown graphically in Figure 16, which shows that these sequences are preferentially shuttled into exosomes and microvesicles compared to the cells.

Gene Symbol	Transcript ID	Length	Type of RNA	CTX0E0307E H cells	CTX0E0307E H EXO	CTX0E0307E H MV	CTX0E0307E I cells	CTX0E0307E I EXO	CTX0E0307E I MV
AC079949.1	AC079949.1-201	57	Novel miRNA	2629	27006	14873	2425	11433	848
AP000318.1	AP000318.1-201	64	Novel miRNA	1353	9379	11002	7469	2963	419
AL161626.1	AL161626.1-201	57	Novel miRNA	471	4450	3712	291	1263	129
AC004943.1	AC004943.1-201	81	Novel miRNA	24	81	43	23	94	5
AL121897.1	AL121897.1-201	89	Novel miRNA	6	22	14	2	30	3

Table 12: Identification of putative novel miRNA sequences using GENCODE in exosomes (EXO), microvesicles (MV) and producer cells. CTX0E03 07EI MV reads are misrepresented due to the lower amount of starting material (table 11). The transcript IDs are taken from the Ensembl database (www.ensembl.org).

Validation and of novel miRNAs

[0298] AC079949.1-201 (SEQ ID NO:738)

Gene: AC079949.1 ENSG00000239776

>12 dna:chromosome chromosome:GRCh37:12:127650616:127650672:1

GGCCGCGCCCCGTTTCCAGGACAAAGGGCACTCCGCACCGGACCCTGGTCCCAGCG

[0299] For AC079949.1-201 putative mature miRNA, **gaccaggguccggugcgagug (SEQ ID NO:745)** was identified as the possible 5' stem mature miRNA using <http://mirna.imbb.forth.gr/MatureBayes.html>, a tool for finding mature miRNA within a miRNA precursor sequence using a Naive Bays classifier. Its presence validation was performed using AG-GGTCCGGTGCGGAGT (SEQ ID NO:746) primer sequence. This sequence was entered in mirbase (<http://www.mirbase.org/>) and the following miRNA was found with similar sequence: *Bos taurus* miR-2887-1 (Accession No. MIMAT0013845).

[0300] bta-miR-2887 : 9-20 (SEQ ID NO:747)

AC079949 (5)	2	ggguccggugcg	13
bta-miR-2887	9	ggguccggugcg	20

[0301] The presence of this novel miRNA was tested by qRT-PCR on purified exosomes retro transcribed miRNA.

[0302] The same analysis was performed using the 3' stem of AC079949, sequence TGCGGAGTGCCCTTTGTCTT (SEQ ID NO:748), but in this case no similar miRNA was identified in mirbase.

[0303] AP000318.1-201 (SEQ ID NO:739)

Gene: AP000318.1 ENSG00000266007

>21 dna:chromosome chromosome:GRCh37:21:35677430:35677493:1

CCCACTCCCTGGCGCCGCTTGTGGAGGGCCCAAGTCCTTCTGATTGAGGCCCAACCCGTGGAAG

[0304] For AP000318.1-201 putative mature miRNA, **ggagggcccaaguccuucugau (SEQ ID NO:744)** was identified as the possible 5' stem mature miRNA. Its presence validation was performed using GGAGGGCCCAAGTCCTTCTGAT (SEQ ID NO:749) primer sequence. *Caenorhabditis remanei* miR-55 stem-loop was identified as similar miRNA. Primer validation was again carried out by qRT-PCR.

crm-miR-55-5p : 4-17 (SEQ ID NO:750)

EP 3 033 418 B1

AP000318.1 20 cccaaguccuucug 7
 ||||| |||||
 crm-miR-55-5p 4 cccaagugcuucug 17

[0305] AL161626.1-201 (SEQ ID NO:740)

Gene: AL161626.1 ENSG00000241781

>9 dna:chromosome chromosome:GRCh37:9:79186731:79186787:1

CGCCGGGACCGGGGTCCGGGGCGGAGTGCCCTTCCTCCTGGGAAACGGGGTGCGGC

[0306] For AL161626.1-201 putative mature miRNA, **ggcggagugcccuuccug** (SEQ ID NO:743) was identified as the possible 5' stem mature miRNA. Its presence validation was performed using CGGAGTGCCCTTCCTCCT (SEQ ID NO:751) primer sequence. *Zea mays* miR164c stem-loop and *Achypodium distachyon* miR164f stem-loop were identified as similar miRNA. Primer validation was again carried out by qRT-PCR.

zma-miR164c-3p : 4-15 (SEQ ID NO:752)

AL161626.1 5 gugcccuucuc 16
 |||||
 zma-miR164c-3p 4 gugcccuucuc 15

[0307] AC004943.1 (SEQ ID NO:741)

Gene: AC004943.1 ENSG00000265573

>16 dna:chromosome chromosome:GRCh37:16:72821592:72821672:-1

GCTTCACGTCCCCACCGCGGCGGCGGCGGTGGCAGTGGCGGCGGCGGCGGCGGTGGCGG
 CGGCGGCGGCGGCGGCGGCGCTC

[0308] AL121897.1 (SEQ ID NO:742)

Gene: AL121897.1 ENSG00000264308

>20 dna:chromosome chromosome:GRCh37:20:30865503:30865591:1

GCCGCCCCCGCCGCCGCCGCCGCCGCCGCCGCCGCCGCCGCCGCCGCTTTCGGCTCGGG
 CCTCAGGTGAGTCGGAGGGGCCGGGCGCC

Miscellaneous RNA (misc_RNA), including novel putative

[0309] Misc_RNA is short for miscellaneous RNA, a general term for a series of miscellaneous small RNA. Miscellaneous transcript feature are not defined by other RNA keys.

[0310] List of top ranking previously known and novel misc_RNAs identified using GENCODE sequence data set:

Gene Symbol	Transcript ID	Length	Type of RNA	CTX0E0307EH	CTX0E0307EH	CTX0E0307EH	CTX0E0307EH	CTX0E0307EH	CTX0E0307EH
				cells	EXO	MV	cells	EXO	MV
RPPH1	RPPH1-201	333	misc RNA	76	2229	1785	0	1077	197
RMRP	RMRP-201	264	misc RNA	139	1803	1443	191	659	87
RPPH1	RPPH1-001	638	misc RNA	182	931	1372	795	2017	157
VTRNA1-1	VTRNA1-1-201	99	misc RNA	43	720	52	247	210	9
Y_RNA	Y_RNA.321-201	93	Novel misc RN	159	196	661	960	903	217
Y_RNA	Y_RNA.725-201	95	Novel misc RN	1092	18	74	1005	39	11
Y_RNA	Y_RNA.125-201	96	Novel misc RN	1079	15	58	906	27	12
Y_RNA	Y_RNA.118-201	99	Novel misc RN	134	12	9	156	45	7
Y_RNA	Y_RNA.394-201	109	Novel misc RN	9	9	7	33	13	1
Y_RNA	Y_RNA.687-201	111	Novel misc RN	36	6	15	103	41	10
Y_RNA	Y_RNA.144-201	102	Novel misc RN	129	5	21	187	84	5
Y_RNA	Y_RNA.337-201	105	Novel misc RN	7	4	0	15	4	0
Y_RNA	Y_RNA.413-201	97	Novel misc RN	136	4	8	125	46	3
Y_RNA	Y_RNA.30-201	103	Novel misc RN	74	3	3	62	21	2

Table 13: Identification of misc_RNA, including putative novel misc_RNA, sequences using GENCODE in exosomes (EXO), microvesicles (MV) and producer cells. (CTX0E03 07EI MV reads are misrepresented due to the lower amount of starting material - Table 11). The transcript IDs are taken from the Ensembl database (www.ensembl.org).

[0311] Among the misc_RNA the following sequences were found preferentially down or up shuttled in exosomes and MV: RPHI, RMRP, and VTRNA1-1 up shuttled and Y_RNA.725-201, and Y_RNA.125-201 down respectively. RPHI is a ribonuclease P RNA component H1. RMRP gene encodes the RNA component of mitochondrial RNA processing endoribonuclease, which cleaves mitochondrial RNA at a priming site of mitochondrial DNA replication. This RNA also interacts with the telomerase reverse transcriptase catalytic subunit to form a distinct ribonucleoprotein complex that has RNA-dependent RNA polymerase activity and produces double-stranded RNAs that can be processed into small interfering RNA. VTRNA1-1 is vault RNA component 1. Vaults are large cytoplasmic ribonucleoproteins and they are composed of a major vault protein, MVP, 2 minor vault proteins, TEP1 and PARP4, and a non-translated RNA component, VTRNA1-1. Y_RNA.725-201, and Y_RNA.125-201 are novel misc_RNAs and their function is not defined.

Metazoa miscellaneous RNA

[0312] The signal recognition particle RNA, also known as 7SL, 6S, ffs, or 4.5S RNA, is the RNA component of the signal recognition particle (SRP) ribonucleoprotein complex. SRP is a universally conserved ribonucleoprotein that directs the traffic of proteins within the cell and allows them to be secreted. The SRP RNA, together with one or more SRP proteins contributes to the binding and release of the signal peptide. The RNA and protein components of this complex are highly conserved but do vary between the different kingdoms of life.

[0313] List of top ranking Metazoa misc_RNAs identified using GENCODE sequence data set:

Gene Symbol	Transcript ID	Length	Type of RNA	CTX0E0307EH	CTX0E0307EH	CTX0E0307EH	CTX0E0307EH	CTX0E0307EH	CTX0E0307EH
				cells	EXO	MV	cells	EXO	MV
Metazoa_SRP	Metazoa_SRP.791-201	288	Metazoan signal recogn	679	2324	2058	771	2698	465
Metazoa_SRP	Metazoa_SRP.561-201	294	Metazoan signal recogn	634	2006	1683	744	2147	432
Metazoa_SRP	Metazoa_SRP.864-201	297	Metazoan signal recogn	252	1884	1544	78	170	148
Metazoa_SRP	Metazoa_SRP.824-201	297	Metazoan signal recogn	438	881	958	505	1860	342
Metazoa_SRP	Metazoa_SRP.72-201	278	Metazoan signal recogn	441	630	631	494	2184	349
Metazoa_SRP	Metazoa_SRP.151-201	307	Metazoan signal recogn	377	464	470	432	1431	265
Metazoa_SRP	Metazoa_SRP.208-201	277	Metazoan signal recogn	382	410	431	422	1104	242
Metazoa_SRP	Metazoa_SRP.501-201	280	Metazoan signal recogn	265	272	266	236	434	44
Metazoa_SRP	Metazoa_SRP.682-201	298	Metazoan signal recogn	12	52	21	10	13	2

Table 14: Identification signal recognition particle RNA (misc_RNA) sequences using GENCODE in exosomes (EXO), microvesicles (MV) and producer cells. The transcript IDs are taken from the Ensembl database (www.ensembl.org).

RRNA (ribosomal RNA)

[0314] Ribosomal RNA (rRNA) forms part of the protein-synthesizing organelle known as a ribosome and that is exported to the cytoplasm to help translate the information in messenger RNA (mRNA) into protein. Eukaryotic ribosome (80S) rRNA components are: large unit (rRNA 5S, 5.8S, and 28S) small unit (rRNA 18S). Both rRNA 28S and 5.8S are selectively up-shuttled in exosomes and MV.

[0315] List of top ranking rRNA identified using GENCODE sequence data set:

Gene Symbol	Transcript ID	Length	Type of RNA	CTX0E0307EH cells	CTX0E0307EH EXO	CTX0E0307EH MV	CTX0E0307EI cells	CTX0E0307EI EXO	CTX0E0307EI MV
RNA5-8SP6	RNA5-8SP6-201	152	rRNA	205008	1148190	706558	213187	135909	14732
RNA28S5	RNA28S5-001	432	rRNA	86111	458585	516754	62829	390237	47483
RNA18S5	RNA18S5-001	599	rRNA	74634	52055	61639	116874	138484	14616
RNA5-8SP2	RNA5-8SP2-201	152	rRNA	6488	1719	1540	9231	3112	149
RNA5-8SP5	RNA5-8SP5-201	152	rRNA	2794	7393	3924	7314	3579	232

Table 15: Identification rRNA sequences using GENCODE in exosomes (EXO), microvesicles (MV) and producer cells. The transcript IDs are taken from the Ensembl database (www.ensembl.org).

Small nucleolar RNA: snoRNA

[0316] Small nucleolar RNAs (snoRNAs) are a class of small RNA molecules that primarily guides chemical modifications of other RNAs, mainly ribosomal RNAs, transfer RNAs and small nuclear RNAs. There are two main classes of snoRNA, the C/D box snoRNAs which are associated with methylation, and the H/ACA box snoRNAs which are associated with pseudouridylation.

[0317] List of top ranking snoRNA identified using GENCODE sequence data set:

Gene Symbol	Transcript ID	Length	Type of RNA	CTX0E0307EH cells	CTX0E0307EH EXO	CTX0E0307EH MV	CTX0E0307EI cells	CTX0E0307EI EXO	CTX0E0307EI MV
SNORD3A	SNORD3A-201	216	snoRNA	1433	2085	1621	906	1732	120
SNORD3C	SNORD3C-201	216	snoRNA	1169	1702	1220	639	1176	86
SNORD29	SNORD29-201	65	snoRNA	28130	1633	1070	36677	1752	45
SNORD83B	SNORD83B-201	93	snoRNA	1835	675	487	638	575	29
SNORD30	SNORD30-201	70	snoRNA	29743	254	244	29071	283	24

Table 16: Identification of snoRNA sequences using GENCODE in exosomes (EXO), microvesicles (MV) and producer cells. The transcript IDs are taken from the Ensembl database (www.ensembl.org).

Small nuclear RNA (snRNA)

[0318] Small nuclear ribonucleic acid (snRNA), also commonly referred to as U-RNA, is a class of small RNA molecules that make up the major spliceosome are named U1, U2, U4, U5, and U6, and participate in several RNA-RNA and RNA-protein interactions. Their primary function is in the processing of pre-mRNA (hnRNA) in the nucleus. They have also been shown to aide in the regulation of transcription factors (7SK RNA) or RNA polymerase II (B2 RNA), and maintaining the telomeres.

[0319] List of top ranking snRNA identified using GENCODE sequence data set:

Gene Symbol	Transcript ID	Length	Type of RNA	CTX0E0307EH cells	CTX0E0307EH EXO	CTX0E0307EH MV	CTX0E0307EI cells	CTX0E0307EI EXO	CTX0E0307EI MV
U2	U2.38-201	191	snRNA	1354	71596	49223	751	35290	1919
U2	U2.6-201	192	snRNA	834	15561	13594	303	8146	272
U1	U1.81-201	164	snRNA	584	10901	7307	91	3197	121
U1	U1.90-201	167	snRNA	533	9927	6689	48	2187	84
U2	U2.7-201	191	snRNA	201	9267	3109	288	6736	262

Table 17A: Identification of snRNA sequences using GENCODE in exosomes (EXO), microvesicles (MV) and producer cells. The transcript IDs are taken from the Ensembl database (www.ensembl.org).

LincRNA and novel lincRNA

[0320] Large intergenic non-coding RNAs (lincRNAs) are emerging as key regulators of diverse cellular processes. Determining the function of individual lincRNAs remains a challenge. Long non-coding RNAs (long ncRNAs, lincRNA) are non-protein coding transcripts longer than 200 nucleotides.

[0321] List of top ranking previously known and novel lincRNAs identified using GENCODE sequence data set:

Gene Symbol	Transcript ID	Length	Type of RNA	CTX0E0307EH cells	CTX0E0307EH EXO	CTX0E0307EH MV	CTX0E0307EI cells	CTX0E0307EI EXO	CTX0E0307EI MV
RP11-108M9.3	RP11-108M9.3-001	1761	Novel lincRNA	244	159	240	539	324	45
RP11-329L6.1	RP11-329L6.1-001	507	Novel lincRNA	19	70	41	29	84	2
RP11-160E2.6	RP11-160E2.6-001	637	Novel lincRNA	228	67	115	489	74	6
AC004528.3	AC004528.3-001	107	Novel lincRNA	16	58	46	14	55	4
MALAT1	MALAT1-201	4585	lincRNA	150	308	234	26	182	12
GAS5	GAS5-007	2743	lincRNA	12024	215	120	46501	875	13

Table 17B: Identification of lincRNA and putative novel lincRNA sequences using GENCODE in exosomes (EXO), microvesicles (MV) and producer cells. The transcript IDs are taken from the Ensembl database (www.ensembl.org).

[0322] GAS5 lincRNA is highly expressed in cell producer compared to in exosomes and microvesicles (down shuttled in both exosomes and MV).

mRNA

[0323] Coding sequencing mRNA were also identified.

Gene Symbol	Transcript ID	Length	Type of RNA	CTX0E0307EH cells	CTX0E0307EH EXO	CTX0E0307EH MV	CTX0E0307EI cells	CTX0E0307EI EXO	CTX0E0307EI MV
EEF2	EEF2-201	9407	mRNA	710	578	449	1155	471	33
MTRNR2L8	MTRNR2L8-201	1290	mRNA	1383	548	642	1323	258	15
NES	NES-001	8635	mRNA	668	406	234	1448	267	20
VIM	VIM-001	8316	mRNA	563	911	501	1500	618	36

Table 18: Identification of mRNA sequences using GENCODE in exosomes (EXO), microvesicles (MV) and producer cells. The transcript IDs are taken from the Ensembl database (www.ensembl.org).

Example 17A-C: Conclusion

[0324] The main scope of the deep sequence analysis was to identify their miRNA components in neural stem cell-derived vesicles (exosomes and microvesicles). This analysis identified a new set of known and novel miRNAs that are

preferentially shuttled into both exosomes and MV. Among the identified miRNAs already included in mirbase database were hsa-miR-1246, hsa-miR-4488, hsa-miR-4492, hsa-miR-4508, hsa-miR-4516, hsa-miR-4532, and among the novel miRNAs were AC079949.1, AP000318.1, AL161626.1, AC004943.1, AL121897.1. Top ranking shuttled miRNAs, including novel ones were validated by qRT-PCR in exosomes.

[0325] The size distribution of shuttle RNA, as shown here, is mostly in the range of 20 to 200 nt and other RNA species are released by cells into the extracellular space. By deep sequencing and GENCODE sequence set analysis we found a greater complexity and diversity of non-coding RNA transcripts. We extended this analysis with detailed evaluation and this led to the discovery of preferentially up (defined as log2 fold change ≥ 2) and down (defined as log2 fold change ≤ -2) shuttle of other non-coding RNAs in both exosomes and microvesicles. Differentially shuttled non coding RNA were found in almost all the non-coding RNA subtypes, ribosomal RNA (rRNA), small nucleolar (snoRNA), small nuclear RNA (snRNA), microRNA (miRNA), miscellaneous other RNA (misc_RNA, e.g. RMRP, vault RNA, metazoa SRP, and RNY), and large intergenic non-coding RNAs (lincRNAs).

[0326] The unequal distribution of the detected RNA species over cellular and shuttle RNA, combined with increasing evidence for their role in gene regulation strongly suggest that cells specifically release these RNAs to modify the function of target cells.

D) Deep sequencing of CTX0E03 cell and exosome miRNA expression from 6 week bioreactor culture

[0327] Next generation deep sequencing was also carried out on CTX0E03 cells and their derived exosomes, following culture for six weeks in an Integra bioreactor. The results showed that hsa-miR-1246, hsa-miR-4492, hsa-miR-4532, and hsa-miR-4488 are also up-shuttled in exosomes derived from 6 week Integra culture (6W). In EXO 6W a total of 61 miRNA types are up-shuttled. Up-shuttled miRNAs with more than 250 reads per exosome sample are listed in Figure 13E.

Conclusions

[0328] Hsa-miR-1246, hsa-miR-4492, hsa-miR-4532, and hsa-miR-4488 are still up-shuttled in EXO 6W as observed on proliferative EXO (07EI & EH; Figure 13 A&B). New up-shuttled miRNAs are also identified, including hsa-miR-4792.

[0329] 20.53 % of the identified miRNA are up-shuttled in the exosomes derived from 6 week Integra CTX cultures (shown in Figure 13C, middle panel). This compares to 99% of the identified miRNAs that are up-shuttled in the exosomes derived from proliferative CTX0E03 cultures (Figure 13C, top panel).

E) Deep sequencing of CTX0E03 cell and exosome miRNA expression from 11 week bioreactor culture

[0330] Next generation deep sequencing was also carried out on CTX0E03 cells and their derived exosomes, following culture for 11 weeks in an Integra bioreactor. Three samples were tested.

[0331] In sample 1, 9 miRNA species are up-shuttled, all of which have more than 250 reads, as shown in Figure 13F.

[0332] In sample 2, 68 miRNA species are up-shuttled into the exosomes. The miRNAs with more than 250 reads per exosome sample are shown in Figure 13G.

[0333] In sample 3, 47 miRNA species are up-shuttled. Figure 13H shows the three miRNA species with a read count >250: hsa-miR-10b-5p, hsa-miR-1246 and hsa-miR-486-5p.

Conclusions

[0334]

Table E1: W11 summary table of reads and log2 values of miRNA types previously reported as up-shuttled in proliferative CTX0E03 exosomes.

MiRNA	Reads						Log2		
	Cell1	Cell2	Cell3	EXO1	EXO2	EXO3	EXO1	EXO2	EXO3
hsa-miR-4488	0	0	1	0	0	0	N/A	N/A	N/A
hsa-miR-4492	0	1	0	1	0	0	N/A	N/A	N/A

EP 3 033 418 B1

(continued)

MiRNA	Reads						Log2		
	Cell1	Cell2	Cell3	EXO1	EXO2	EXO3	EXO1	EXO2	EXO3
hsa-miR-4532	0	0	0	0	0	0	N/A	N/A	N/A
hsa-miR-1246	483	1122	3470	18	2726	24152	-4.20	1.26	2.99

[0335] Hsa-miR-1246 is present in 11W exosomes, but was only observed to be up-shuttled in EXO3.

[0336] Hsa-miR-4488, hsa-miR-4492, and hsa-miR-4532, identified in proliferative CTX0E03 cells and their exosomes, are almost absent in 11 week samples (both cells and exosomes).

[0337] Hsa-miR-486-5p was the only miRNA up-shuttled in all three EXO W11 samples.

[0338] An average 12.22% of the identified miRNAs are up-shuttled in the exosomes derived from 11 week Integra CTX0E03 cultures (Figure 13C, lower panel).

Comparative summary tables

Comparative analysis of miRNA expression in EXO samples sorted by largest reads in EXO derived from proliferative CTX0E03/07EH

[0339]

MiRNA	Reads						LOG2			
	EXO 07EH	EXO 07EI	EXO WK6	EXO1 WK11	EXO2 WK11	EXO3 WK11	EXO WK6	EXO1 WK11	EXO2 WK11	EXO3 WK11
hsa-miR-1246	111092	83958	77678	18	2726	24152	-1.40	-13.43	-9.32	-5.96
hsa-miR-4492	5188	22482	9528	1			-2.65	-15.82		
hsa-miR-4532	3368	6419	3463				-2.39			
hsa-miR-4488	1389	20618	11048				-2.12			
hsa-miR-4485	386	43	13		11	69	-5.46		-8.80	-5.93
hsa-miR-4508	188	723	1245			20	-0.68			-9.51
hsa-miR-4516	135	904	1153				-0.96			
hsa-miR-4497	73	28	61				-1.64			
hsa-miR-1973	50	5	1				-7.05			
hsa-miR-3195	48	17	14				-3.40			
hsa-miR-4466	43	12	15			9	-2.38			-5.99
hsa-let-7a-5p	20	15	18306	69521	496463	509272	7.80	9.77	9.47	9.72
hsa-miR-99b-5p	19	15	25938	6984	56360	50611	8.35	6.50	6.37	6.44
hsa-miR-21-5p	18	41	35754	12191	134432	150516	8.01	6.50	6.82	7.21
hsa-miR-92a-3p	18	7	15054	5504	89826	73306	8.00	6.59	7.48	7.41
hsa-miR-3676-5p	17	174	59				-2.88			
hsa-miR-4792	15	24	254				1.54			
hsa-miR-664-5p	13	8	1				-5.53			
hsa-miR-100-5p	11	5	30124	22	56751	43080	9.68	-0.69	7.50	7.32
hsa-miR-1291	10	17	18			34	-1.83			-3.79
hsa-miR-16-5p	10	4	5502	15	6580	9510	7.32	-1.15	4.48	5.24
hsa-miR-4284	10	11	5				-3.26			
hsa-miR-663b	9	26	39				-1.03			
hsa-miR-25-3p	8		1383	5825	42402	29320	6.09	8.21	7.94	7.63
hsa-miR-3656	8	10	743	2023	1		4.20	5.69	-8.43	
hsa-miR-181a-5p	7	1	81465	56064	1289686	1023049	12.19	11.69	13.08	12.96
hsa-miR-3654	6	13	3			9	-3.90			-5.19
hsa-miR-26a-5p	6	17	18836	3542	49957	54929	8.45	6.08	6.76	7.12

Table E2: Summary table listing miRNA reads and log2. Log2 is calculated using the normalized ratio of either EXO 6W or EXO 11W samples / averaged reads in EXO derived from proliferative cells. Up-shuttled miRNAs ($\log_2 > 2$), in EXO derived from CTX0E03 cultured for 6 and 11 weeks in Integra flasks, are highlighted in lighter grey and down-shuttled ($\log_2 < 2$) in darker grey respectively. The table presents only the top 30 more abundant miRNAs.

Comparative analysis of miRNA expression in EXO samples sorted by largest reads in EXO derived from 6W Integra CTX0E03 culture

[0340]

EP 3 033 418 B1

MiRNA	Reads						LOG2			
	EXO 07EH	EXO 07EI	EXO WK6	EXO1 WK11	EXO2 WK11	EXO3 WK11	EXO WK6	EXO1 WK11	EXO2 WK11	EXO3 WK11
hsa-miR-181a-5p	7	1	81465	56064	1289686	1023049	12.19	11.69	13.08	12.96
hsa-miR-1246	111092	83958	77678	18	2726	24152	-1.40	-13.43	-9.32	-5.96
hsa-miR-127-3p		6	36828	29	22175	48371	11.39	1.13	7.57	8.92
hsa-miR-21-5p	18	41	35754	12191	134432	150516	8.01	6.50	6.82	7.21
hsa-miR-100-5p	11	5	30124	22	56751	43080	9.68	-0.69	7.50	7.32
hsa-miR-99b-5p	19	15	25938	6984	56360	50611	8.35	6.50	6.37	6.44
hsa-miR-26a-5p	6	17	18836	3542	49957	54929	8.45	6.08	6.76	7.12
hsa-let-7a-5p	20	15	18306	69521	496463	509272	7.80	9.77	9.47	9.72
hsa-miR-191-5p	4	1	17845	10674	130737	104349	10.49	9.79	10.27	10.17
hsa-miR-27b-3p	1		15232	7344	33638	47061	12.55	11.55	10.60	11.31
hsa-miR-92a-3p	18	7	15054	5504	89826	73306	8.00	6.59	7.48	7.41
hsa-miR-146b-5p		1	12960	12	29712	13822	12.47	2.44	10.58	9.69
hsa-let-7i-5p	1	2	12914	6928	157331	149116	10.83	9.98	11.34	11.49
hsa-miR-222-3p	3	2	12671	46	89141	64302	10.03	1.97	9.75	9.50
hsa-let-7c	1	2	12414				10.77			
hsa-miR-4488	1389	20618	11048				-2.12			
hsa-miR-4492	5188	22482	9528	1			-2.65	-15.82		
hsa-let-7b-5p	5	6	8914	34	86886	33259	8.40	0.41	8.59	7.43
hsa-let-7f-5p	5	6	8766	13898	293116	427140	8.38	9.09	10.35	11.11
hsa-miR-92b-3p	3	7	8141	73468	420313	291688	8.43	11.65	11.03	10.72
hsa-miR-181b-5p	1	3	7919	6713	226010	167200	9.72	9.53	11.46	11.25
hsa-miR-221-3p	3	3	6391	4927	17702	12768	8.79	8.46	7.17	6.92
hsa-miR-125b-5p	4	4	6107	2181	54267	37113	8.31	6.87	8.37	8.04
hsa-miR-151a-3p	3	4	6014	2776	26457	29135	8.49	7.42	7.53	7.89
hsa-miR-30a-5p	2		5613	15559	136508	96950	10.11	11.63	11.62	11.35
hsa-miR-16-5p	10	4	5502	15	6580	9510	7.32	-1.15	4.48	5.24
hsa-miR-22-3p	4	1	5489	2914	8193	10780	8.79	7.92	6.27	6.89
hsa-miR-99a-5p	1		5116	14	27513	17532	10.98	2.51	10.31	9.88

Table E3: Summary table listing miRNA reads and log2. Log2 is calculated using the normalized ratio of either EXO 6W or EXO 11W samples / averaged reads in EXO derived from proliferative cells. Up-shuttled miRNAs (log2 > 2), in EXO derived from CTX0E03 cultured for 6 and 11 weeks in Integra flasks, are highlighted in lighter grey and down-shuttled (log2 <2) in darker grey respectively. The table presents only the top 30 more abundant miRNAs.

Comparative analysis of miRNA expression in EXO samples sorted by largest reads in EXO3 derived from 11W Integra CTX0E03 culture

[0341]

MiRNA	Reads						LOG2			
	EXO 07EH	EXO 07EI	EXO WK6	EXO1 WK11	EXO2 WK11	EXO3 WK11	EXO WK6	EXO1 WK11	EXO2 WK11	EXO3 WK11
hsa-miR-181a-5p	7	1	81465	56064	1289686	1023049	12.19	11.69	13.08	12.96
hsa-let-7a-5p	20	15	18306	69521	496463	509272	7.80	9.77	9.47	9.72
hsa-let-7f-5p	5	6	8766	13898	293116	427140	8.38	9.09	10.35	11.11
hsa-miR-92b-3p	3	7	8141	73468	420313	291688	8.43	11.65	11.03	10.72
hsa-miR-9-5p	3	3	2196	155707	481660	263043	7.25	13.44	11.93	11.28
hsa-miR-181b-5p	1	3	7919	6713	226010	167200	9.72	9.53	11.46	11.25
hsa-miR-21-5p	18	41	35754	12191	134432	150516	8.01	6.50	6.82	7.21
hsa-let-7i-5p	1	2	12914	6928	157331	149116	10.83	9.98	11.34	11.49
hsa-miR-191-5p	4	1	17845	10674	130737	104349	10.49	9.79	10.27	10.17
hsa-miR-30a-5p	2		5613	15559	136508	96950	10.11	11.63	11.62	11.35
hsa-miR-92a-3p	18	7	15054	5504	89826	73306	8.00	6.59	7.48	7.41
hsa-miR-10a-5p		1	3754	5840	132775	66316	10.69	11.37	12.74	11.96
hsa-miR-222-3p	3	2	12671	46	89141	64302	10.03	1.97	9.75	9.50
hsa-let-7c-5p				13106	68120	59498				
hsa-miR-26a-5p	6	17	18836	3542	49957	54929	8.45	6.08	6.76	7.12
hsa-miR-423-5p	1	1	538	2606	45587	52945	6.81	9.13	10.12	10.55
hsa-miR-99b-5p	19	15	25938	6984	56360	50611	8.35	6.50	6.37	6.44
hsa-miR-127-3p		6	36828	29	22175	48371	11.39	1.13	7.57	8.92
hsa-miR-27b-3p	1		15232	7344	33638	47061	12.55	11.55	10.60	11.31
hsa-miR-100-5p	11	5	30124	22	56751	43080	9.68	-0.69	7.50	7.32
hsa-let-7e-5p		2	1932	5398	27921	41799	8.73	10.26	9.49	10.29
hsa-miR-486-5p	3	1	3503	44199	20054	37232	8.47	12.17	7.89	9.01
hsa-miR-125b-5p	4	4	6107	2181	54267	37113	8.31	6.87	8.37	8.04
hsa-let-7b-5p	5	6	8914	34	86886	33259	8.40	0.41	8.59	7.43
hsa-miR-182-5p	3	15	833	22	33064	32491	4.32	-0.88	6.53	6.73
hsa-miR-30d-5p		2	3946	22	26066	32148	9.76	2.32	9.39	9.91
hsa-miR-25-3p	8		1383	5825	42402	29320	6.09	8.21	7.94	7.63
hsa-miR-151a-3p	3	4	6014	2776	26457	29135	8.49	7.42	7.53	7.89

Table E4: Summary table listing miRNA reads and log2. Log2 is calculated using the normalized ratio of either EXO 6W or EXO 11W samples / averaged reads in EXO derived from proliferative cells. Up-shuttled miRNAs ($\log_2 > 2$), in EXO derived from CTX0E03 cultured for 6 and 11 weeks in Integra flasks, are high-lighted in lighter grey and down-shuttled ($\log_2 < 2$) in darker grey respectively. The table presents only the top 30 more abundant miRNAs.

Conclusions for comparative summary of miRNA reads present in exosome samples

[0342] Hsa-miR-1246, hsa-miR-4492, hsa-miR-4532, and hsa-miR-4488 are the most up-shuttled miRNA types in exosomes derived from proliferative CTX0E03 cells.

[0343] Hsa-miR-1246, hsa-miR-4492, hsa-miR-4532, and hsa-miR-4488 are still present in EXO 6W sample, but hsa-miR-4492, hsa-miR-4532, and hsa-miR-4488 are almost absent in EXO 11W samples.

[0344] Hsa-miR-181a-5p, hsa-miR-1246, hsa-miR-127-3p, hsa-miR-21-5p, and hsa-miR-100-5p are the top 5 miRNAs present in EXO 6W sample.

[0345] Hsa-miR-181a-5p, hsa-let-7a-5p, hsa-let-7f-5p, hsa-miR-92b-3p, and hsa-miR-9-5p are the top 5 miRNAs present in EXO 11W samples.

Comparative analysis of miRNA expression in cell samples sorted by largest reads in proliferative cell 07EH

[0346]

EP 3 033 418 B1

miRNA	Reads						LOG2			
	Cell (07EH)	Cells (07EI)	Cell 6WK	Cell1 WK11	Cell2 WK11	Cell3 WK11	Cell 6WK	Cell1 WK11	Cell2 WK11	Cell3 WK11
hsa-let-7a-5p	75110	305060	357507	82000	572638	945050	0.08	-0.68	-0.44	-0.58
hsa-miR-92b-3p	9794	242715	51788	7070	95904	269851	-1.46	-2.97	-1.78	-1.15
hsa-miR-21-5p	11943	154626	109713	36935	165335	350731	0.06	-0.15	-0.55	-0.33
hsa-miR-92a-3p	14359	137412	136842	15296	99662	238214	0.41	-1.39	-1.25	-0.86
hsa-miR-127-3p	7064	110806	54214	4015	20215	70609	-0.40	-2.79	-3.03	-2.09
hsa-miR-100-5p	52451	109290	121101	8268	35093	40100	-0.61	-3.12	-3.60	-4.28
hsa-miR-27b-3p	16900	91902	55177	15130	73464	113623	-0.66	-1.16	-1.45	-1.69
hsa-miR-191-5p	12591	89150	58145	16783	74782	151339	-0.37	-0.80	-1.21	-1.06
hsa-miR-26a-5p	9900	88724	164401	24518	127010	271289	1.27	-0.11	-0.30	-0.08
hsa-miR-99b-5p	39457	87399	46207	3177	14278	22554	-1.62	-4.12	-4.52	-4.72
hsa-let-7f-5p	10349	78395	61473	115717	457274	982774	-0.06	2.21	1.63	1.87
hsa-miR-181a-5p	6956	47686	180839	150467	1304081	1909336	2.15	3.25	3.80	3.48
hsa-miR-486-5p	20310	41639	4938	465	1500	4379	-3.85	-5.90	-6.78	-6.10
hsa-miR-30a-5p	2001	35465	16099	23958	94346	201534	-0.47	1.47	0.88	1.11
hsa-miR-98	11760	30440	14559				-1.61			
hsa-miR-151a-3p	2681	29047	15661	5186	22917	58870	-0.42	-0.65	-1.07	-0.57
hsa-miR-21-3p	4089	27733	10626	2682	10247	30683	-1.16	-1.78	-2.42	-1.70
hsa-miR-30d-5p	1977	27307	17522	4802	35022	65318	-0.06	-0.56	-0.26	-0.23
hsa-let-7c	5103	27224	70342				1.43			
hsa-miR-10a-5p	52927	26908	28817	11983	49442	69645	-2.35	-2.25	-2.77	-3.14
hsa-miR-22-3p	1826	26456	10088	5318	33593	54768	-0.79	-0.35	-0.26	-0.42
hsa-miR-182-5p	5531	25885	12376	2528	18519	49416	-1.10	-2.03	-1.72	-1.17
hsa-miR-222-3p	1422	22187	10960	13094	53385	90511	-0.39	1.23	0.69	0.59
hsa-miR-125a-5p	1451	20960	9101	1406	12194	16095	-0.61	-1.94	-1.39	-1.85
hsa-miR-16-5p	2173	19856	69449	8542	42494	81960	2.20	0.54	0.29	0.37
hsa-let-7b-5p	2435	19774	61291	13225	75017	99021	1.96	1.11	1.05	0.58
hsa-miR-151a-5p	1386	19773	10790	1621	9192	17997	-0.28	-1.65	-1.72	-1.61
hsa-let-7e-5p	2449	19035	14175	7948	40674	84912	-0.12	0.41	0.20	0.39

Table E5: Summary table listing miRNA reads and log2. Log2 is calculated using the normalized ratio of either cell 6W or cell 11W samples / averaged reads in proliferative cells. Up-expressed miRNAs (log2 > 2), in CTX0E03 cultured for 6 and 11 weeks in Integra flasks, are highlighted in lighter grey and down-expressed (log2 <2) in darker grey respectively. The table presents only the top 30 more abundant miRNAs.

Comparative analysis of miRNA expression in cell samples sorted by largest reads in cells cultured for 6 week in integra flasks (Cell 6W)

[0347]

EP 3 033 418 B1

miRNA	Reads						LOG2			
	Cell (07EH)	Cells (07EI)	Cell 6WK	Cell1 WK11	Cell2 WK11	Cell3 WK11	Cell 6WK	Cell1 WK11	Cell2 WK11	Cell3 WK11
hsa-let-7a-5p	75110	305060	357507	82000	572638	945050	0.08	-0.68	-0.44	-0.58
hsa-miR-181a-5p	6956	47686	180839	150467	1304081	1909336	2.15	3.25	3.80	3.48
hsa-miR-26a-5p	9900	88724	164401	24518	127010	271289	1.27	-0.11	-0.30	-0.08
hsa-miR-92a-3p	14359	137412	136842	15296	99662	238214	0.41	-1.39	-1.25	-0.86
hsa-miR-100-5p	52451	109290	121101	8268	35093	40100	-0.61	-3.12	-3.60	-4.28
hsa-miR-21-5p	11943	154626	109713	36935	165335	350731	0.06	-0.15	-0.55	-0.33
hsa-let-7c	5103	27224	70342				1.43			
hsa-miR-16-5p	2173	19856	69449	8542	42494	81960	2.20	0.54	0.29	0.37
hsa-let-7f-5p	10349	78395	61473	115717	457274	982774	-0.06	2.21	1.63	1.87
hsa-let-7b-5p	2435	19774	61291	13225	75017	99021	1.96	1.11	1.05	0.58
hsa-miR-191-5p	12591	89150	58145	16783	74782	151339	-0.37	-0.80	-1.21	-1.06
hsa-miR-27b-3p	16900	91902	55177	15130	73464	113623	-0.66	-1.16	-1.45	-1.69
hsa-miR-127-3p	7064	110806	54214	4015	20215	70609	-0.40	-2.79	-3.03	-2.09
hsa-miR-92b-3p	9794	242715	51788	7070	95904	269851	-1.46	-2.97	-1.78	-1.15
hsa-miR-9-5p	518	7957	49128	113170	696970	1266994	3.25	5.82	5.87	5.87
hsa-miR-99b-5p	39457	87399	46207	3177	14278	22554	-1.62	-4.12	-4.52	-4.72
hsa-miR-146b-5p	4552	8434	38088	7472	43047	20961	1.29	0.30	0.26	-1.64
hsa-miR-125b-5p	1002	17965	34339	2841	21418	24635	1.61	-0.62	-0.28	-0.94
hsa-miR-124e	3973	1783	32042	483	1122	3470	1.56	-3.13	-4.48	-3.72
hsa-miR-10a-5p	52927	26908	28817	11983	49442	69645	-2.35	-2.25	-2.77	-3.14
hsa-let-7i-5p	3015	17802	21469	33910	167628	432785	0.40	2.43	2.16	2.67
hsa-miR-99a-5p	773	2767	19124	1826	11313	14604	2.54	0.51	0.58	0.08
hsa-miR-30d-5p	1977	27307	17522	4802	35022	65318	-0.06	-0.56	-0.26	-0.23
hsa-let-7g-5p	959	15467	16691	4240	19981	36131	0.75	0.14	-0.20	-0.21
hsa-miR-30a-5p	2001	35465	16099	23958	94346	201534	-0.47	1.47	0.88	1.11
hsa-miR-151a-3p	2681	29047	15661	5186	22917	58870	-0.42	-0.65	-1.07	-0.57
hsa-miR-98	11760	30440	14559				-1.61			
hsa-miR-204-5p	113	1378	14322	1554	8147	18192	3.91	2.07	1.89	2.18

Table E6: Summary table listing miRNA reads and log2. Log2 is calculated using the normalized ratio of either cell 6W or cell 11W samples / averaged reads in proliferative cells. Up-expressed miRNAs (log2 > 2), in CTX0E03 cultured for 6 and 11 weeks in Integra flasks, are high-lighted in lighter grey and down-expressed (log2 < 2) in darker grey respectively. The table presents only the top 30 more abundant miRNAs.

Comparative analysis of miRNA expression in cell samples sorted by largest reads in cells cultured 11 week in Integra flasks (Cell1 W11)

[0348]

miRNA	Reads						LOG2			
	Cell (07EH)	Cells (07EI)	Cell 6WK	Cell1 WK11	Cell2 WK11	Cell3 WK11	Cell 6WK	Cell1 WK11	Cell2 WK11	Cell3 WK11
hsa-miR-181a-5p	6956	47686	180839	150467	1304081	1909336	2.15	3.25	3.80	3.48
hsa-miR-9-5p	518	7957	49128	113170	696970	1266994	3.25	5.82	5.87	5.87
hsa-let-7f-5p	10349	78395	61473	115717	457274	982774	-0.06	2.21	1.63	1.87
hsa-let-7a-5p	75110	305060	357507	82000	572638	945050	0.08	-0.68	-0.44	-0.58
hsa-let-7i-5p	3015	17802	21469	33910	167628	432785	0.40	2.43	2.16	2.67
hsa-miR-21-5p	11943	154626	109713	36935	165335	350731	0.06	-0.15	-0.55	-0.33
hsa-miR-181b-5p	1382	12606	11845	17620	135852	275865	0.30	2.24	2.62	2.77
hsa-miR-26a-5p	9900	88724	164401	24518	127010	271289	1.27	-0.11	-0.30	-0.08
hsa-miR-92b-3p	9794	242715	51788	7070	95904	269851	-1.46	-2.97	-1.78	-1.15
hsa-miR-92a-3p	14359	137412	136842	15296	99662	238214	0.41	-1.39	-1.25	-0.86
hsa-miR-30a-5p	2001	35465	16099	23958	94346	201534	-0.47	1.47	0.88	1.11
hsa-miR-191-5p	12591	89150	58145	16783	74782	151339	-0.37	-0.80	-1.21	-1.06
hsa-let-7c-5p				18010	106589	130919				
hsa-miR-27b-3p	16900	91902	55177	15130	73464	113623	-0.66	-1.16	-1.45	-1.69
hsa-let-7b-5p	2435	19774	61291	13225	75017	99021	1.96	1.11	1.05	0.58
hsa-miR-423-5p	1080	8893	3219	9325	40672	93044	-1.13	1.77	1.33	1.65
hsa-miR-222-3p	1422	22187	10960	13094	53385	90511	-0.39	1.23	0.69	0.59
hsa-let-7e-5p	2449	19035	14175	7948	40674	84912	-0.12	0.41	0.20	0.39
hsa-miR-16-5p	2173	19856	69449	8542	42494	81960	2.20	0.54	0.29	0.37
hsa-miR-25-3p	436	5303	8387	5633	31971	77502	1.19	1.98	1.92	2.33
hsa-miR-127-3p	7064	110806	54214	4015	20215	70609	-0.40	-2.79	-3.03	-2.09
hsa-miR-10a-5p	52927	26908	28817	11983	49442	69645	-2.35	-2.25	-2.77	-3.14
hsa-miR-30d-5p	1977	27307	17522	4802	35022	65318	-0.06	-0.56	-0.26	-0.23
hsa-miR-151a-3p	2681	29047	15661	5186	22917	58870	-0.42	-0.65	-1.07	-0.57
hsa-miR-22-3p	1826	26456	10088	5318	33593	54768	-0.79	-0.35	-0.26	-0.42
hsa-miR-182-5p	5531	25885	12376	2528	18519	49416	-1.10	-2.03	-1.72	-1.17
hsa-miR-181a-3p	86	803	1186	2490	21343	44745	0.97	3.40	3.93	4.13
hsa-miR-100-5p	52451	109290	121101	8268	35093	40100	-0.61	-3.12	-3.60	-4.28

Table E7: Summary table listing miRNA reads and log2. Log2 is calculated using the normalized ratio of either cell 6W or cell 11W samples / averaged reads in proliferative cells. Up-expressed miRNAs ($\log_2 > 2$), in CTX0E03 cultured for 6 and 11 weeks in Integra flasks, are high-lighted in lighter grey and down-expressed ($\log_2 < 2$) in darker grey respectively. The table presents only the top 30 more abundant miRNAs.

Conclusions for comparative summary of miRNA reads present in cell samples

[0349] Hsa-let-7a-5p, hsa-miR-92b-3p, hsa-miR-21-5p, hsa-miR-92a-3p, and hsa-miR-127-3p are the top 5 most expressed miRNA types in proliferative CTX0E03 cells.

[0350] Hsa-let-7a-5p, hsa-miR-181a-5p, hsa-miR-26a-5p, hsa-miR-92a-3p, hsa-miR-100-5p are the top 5 most expressed miRNA types in CTX0E03 Integra 6W culture.

[0351] Hsa-miR-181a-5p, hsa-miR-9-5p, hsa-let-7f-5p, hsa-let-7a-5p, and hsa-let-7i-5p are the top 5 most expressed miRNA types in CTX0E03 Integra 11W cultures.

[0352] Hsa-miR-181a-5p and hsa-miR-9-5p are up-expressed in all cell samples cultured in Integra flasks (6 and 11 weeks).

[0353] Hsa-let-7i-5p, hsa-let-7c-5p, hsa-miR-181a-3p and hsa-miR-181b-5p were solely up-expressed in W11 cells.

[0354] Hsa-miR-181 family seems to play an important role in CTX0E03 long term culture and possible differentiation.

Example 18: Proteomic analysis

Methods

[0355] Exosomes and microvesicle fractions were prepared from a CTX0E03 cell Integra culture (week 2), using differential ultracentrifugation. Exosomes and microvesicles were disrupted in modified RIPA buffer (50mM Tris HCl, pH 8.0, 150mM NaCl, 1% SDS, 0.1% Triton X100, 10mM DTT, 1x Complete protease inhibitor (Roche) and 1x PhosStop phosphatase inhibitor (Roche)) and subjected to manual shearing using a 1mL tuberculin syringe and 25 gauge needle. Samples were re-quantitated post disruption using the Qubit fluorometer (Invitrogen). 20µg of each sample was loaded onto a 4-12% SDS-PAGE gel (Novex, Invitrogen). The gel was excised into forty segments per lane and gel slices were

processed using a robot (ProGest, DigiLab) with the following protocol:

- a) wash with 25mM ammonium bicarbonate followed by acetonitrile;
- b) reduce with 10mM dithiothreitol at 60°C followed by alkylation with 50mM iodoacetamide at room temperature;
- c) digest with trypsin (Promega) at 37°C for 4h;
- d) quench with formic acid;
- e) the supernatant was analysed by mass spectrometry directly without further processing.

Mass Spectrometry

[0356] Each gel digest was analysed by nano LC/MS/MS with a Waters NanoAcquity HPLC system interfaced to a ThermoFisher Q Exactive. Peptides were loaded on a trapping column and eluted over a 75µm analytical column at 350nL/min; both columns were packed with Jupiter Proteo resin (Phenomenex). The mass spectrometer was operated in data-dependent mode, with MS and MS/MS performed in the Orbitrap at 70,000 FWHM and 17,500 FWHM resolution, respectively.

Exosomes

[0357] 2572 proteins were identified by Mass spectrometry in exosomes purified by ultracentrifugation. The exosomes were isolated from the initial stages of an Integra culture (week 2). The gene names and corresponding SWISSPROT accession numbers (in brackets) of all 2572 proteins are listed in Table 19 (in alphabetical order of gene name) and the 100 most abundant proteins are listed in Table 20, in order of decreasing abundance. The characteristic exosome markers CD9, CD81 and Alix (also known as PDCD6IP) are present in the most abundant 100 proteins.

Table 19: Gene names and SWISSPROT accession numbers of all 2572 proteins identified in CTX0E03 exosomes (listed in alphabetical order of gene name).

A1BG (P04217), A2M (P01023), AACS (Q86V21), AAMP (Q13685), AARS (P49588), AARSD1 (Q9BTE6), AASDHPPT (Q9NRN7), ABCA3 (Q99758), ABCE1 (P61221), ABCF1 (Q8NE71), ABCF3 (Q9NUQ8), ABHD10 (Q9NUJ1), ABHD14B (Q96IU4), ABI1 (Q8IZP0), ABR (Q12979), ACAA2 (P42765), ACACA (Q13085), ACADVL (P49748), ACAP2 (Q15057), ACAT1 (P24752), ACAT2 (Q9BWD1), ACBD7 (Q8N6N7), ACLY (P53396), ACO1 (P21399), ACO2 (Q99798), ACOT1 (Q86TX2), ACOT13 (Q9NPJ3), ACOT7 (O00154), ACP1 (P24666), ACSL1 (P33121), ACSL3 (Q95573), ACSL4 (O60488), ACSS2 (Q9NR19), ACTC1 (P68032), ACTG1 (P63261), ACTL6A (O96019), ACTN1 (P12814), ACTN4 (O43707), ACTR10 (Q9NZ32), ACTR1A (P61163), ACTR1 B (P42025), ACTR2 (P61160), ACTR3 (P61158), ADAM10 (O14672), ADAM12 (O43184), ADAMTS15 (Q8TE58), ADAMTS16 (Q8TE57), ADAR (P55265), ADAT2 (Q7Z6V5), ADH5 (P11766), ADI1 (Q9BV57), ADK (P55263), ADRBK1 (P25098), ADRM1 (O16186), ADSL (P30566), ADSS (P30520), AEBP1 (Q8IUX7), AFM (P43652), AGL (P35573), AGRN (O00468), AGT (P01019), AHCY (P23526), AHCYL1 (O43865), AHNK (Q09666), AHSA1 (O95433), AHSG (P02765), AIDA (Q96BJ3), AIFM1 (O95831), AIMP1 (Q12904), AIMP2 (Q13155), AIP (O00170), AK1 (P00568), AK3 (Q9UIJ7), AK4 (P27144), AKAP12 (Q02952), AKAP9 (Q99996), AKR1A1 (P14550), AKR1B1 (P15121), AKR1C1 (Q04828), AKR7A2 (O43488), AKR7A3 (Q95154), AKT1 (P31749), ALCAM (Q13740), ALDH16A1 (Q8IZ83), ALDH3A1 (P30838), ALDH7A1 (P49419), ALDH9A1 (P49189), ALDOA (P04075), ALDOC (P09972), ALKBH2 (Q6NS38), ALKBH4 (Q9NXW9), AMBP (P02760), AMDHD2 (Q9Y303), AMPD2 (Q01433), AMZ2 (Q86W34), ANAPC1 (Q9H1A4), ANAPC4 (Q9UJX5), ANAPC5 (Q9UJX4), ANAPC7 (Q9UJX3), ANKFY1 (Q9P2R3), ANKRD28 (O15084), ANP32A (P39687), ANP32B (Q92688), ANP32E (Q9BTT0), ANXA1 (P04083), ANXA2 (P07355), ANXA4 (P09525), ANXA5 (P08758), ANXA6 (P08133), ANXA7 (P20073), AP1B1 (Q10567), AP1G1 (O43747), AP1M1 (Q9BXS5), AP1S1 (P61966), AP1S2 (P56377), AP2A1 (Q95782), AP2A2 (O94973), AP2B1 (P63010), AP2M1 (Q96CW1), AP2S1 (P53680), AP3B1 (O00203), AP3D1 (O14617), AP3M1 (Q9Y2T2), AP3S1 (Q92572), AP3S2 (P59780), AP4S1 (Q9Y587), APEH (P13798), APEX1 (P27695), API5 (Q9BZZ5), APIP (Q96GX9), APOA1 (P02647), APOA1 BP (Q8NCW5), APOA2 (P02652), APOBEC3C (Q9NRW3), APOC2 (P02655), APOD (P05090),

EP 3 033 418 B1

(continued)

APOH (P02749), APOM (O95445), APPL1 (Q9UKG1), APRT (P07741), AQR (O60306), ARCN1 (P48444), ARF1 (P84077), ARF4 (P18085), ARF5 (P84085), ARF6 (P62330), ARFIP1 (P53367), ARFIP2 (P53365), ARHGAP1 (Q07960), ARHGAP12 (Q8IWW6), ARHGDIA (P52565), ARHGEF1 (Q92888), ARHGEF10 (O15013), ARHGEF7 (Q14155), ARIH1 (Q9Y4X5), ARIH2 (O95376), ARL1 (P40616), ARL2 (P36404), ARL3 (P36405), ARL6IP1 (Q15041), ARL8B (Q9NVJ2), ARMC10 (Q8N2F6), ARMC6 (Q6NXE6), ARMC8 (Q8IUR7), ARMC9 (Q7Z3E5), ARMCX3 (Q9UH62), ARPC1A (Q92747), ARPC1B (O15143), ARPC2 (O15144), ARPC3 (O15145), ARPC4 (P59998), ARPC5 (O15511), ARPC5L (Q9BPX5), ARRDC1 (Q8N5I2), ASB6 (Q9NWX5), ASCC1 (Q8N9N2), ASCC2 (Q9H1I8), ASCC3 (Q8N3C0), ASF1A (Q9Y294), ASH2L (Q9UBL3), ASMTL (O95671), ASNA1 (O43681), ASNS (P08243), ASS1 (P00966), ATG16L1 (Q676U5), ATG3 (Q9NT62), ATG4B (Q9Y4P1), ATG7 (O95352), ATIC (P31939), ATL3 (Q6DD88), ATM (Q13315), ATOX1 (O00244), ATP1A1 (P05023), ATP1B1 (P05026), ATP1B3 (P54709), ATP2B1 (P20020), ATP2B4 (P23634), ATP5B (P06576), ATP5E (P56381), ATP5I (P56385), ATP6AP2 (O75787), ATP6V0D1 (P61421), ATP6V1A (P38606), ATP6V1B2 (P21281), ATP6V1C1 (P21283), ATP6V1D (Q9Y5K8), ATP6V1E1 (P36543), ATP6V1G1 (O75348), ATP6V1 H (Q9UI12), ATR (Q13535), ATRN (O75882), ATXN10 (Q9UBB4), B2M (P61769), B3GAT3 (O94766), B3GNT1 (O43505), B4GALT7 (Q9UBV7), BAG2 (O95816), BAIAP2 (Q9UQB8), BANF1 (O75531), BAT1 (Q13838), BAT3 (P46379), BBOX1 (O75936), BCAS2 (O75934), BCAT1 (P54687), BCCIP (Q9P287), BCL2L13 (Q9BXX5), BCLAF1 (Q9NYF8), BDH2 (Q9BUT1), BICD2 (Q8TD16), BLOC1S1 (P78537), BLVRA (P53004), BLVRB (P30043), BMP1 (P13497), BOLA2 (Q9H3K6), BPGM (P07738), BPHL (Q86WA6), BPNT1 (O95861), BRCC3 (P46736), BRE (Q9NXR7), BROX (Q5VW32), BRP16L (POCB43), BSG (P35613), BST1 (Q10588), BTAF1 (O14981), BUB3 (O43684), BUD31 (P41223), BYSL (Q13895), BZW1 (Q7L1Q6), BZW2 (Q9Y6E2), C10orf119 (Q9BTE3), C10orf58 (Q9BRX8), C10orf76 (Q5T2E6), C11orf54 (Q9H0W9), C11orf68 (Q9H3H3), C12orf10 (Q9HB07), C14orf149 (Q96EM0), C14orf166 (Q9Y224), C15orf58 (Q6ZNW5), C16orf13 (Q96S19), C16orf80 (Q9Y6A4), C1D (Q13901), C1orf123 (Q9NWW4), C1orf50 (Q9BV19), C1orf57 (Q9BSD7), C1RL (Q9NZP8), C20orf11 (Q9NWWU2), C20orf27 (Q9GZN8), C20orf4 (Q9Y312), C21orf59 (P57076), C22orf25 (Q6ICL3), C22orf28 (Q9Y3I0), C2orf29 (Q9UKZ1), C2orf79 (Q6GMV3), C3orf10 (Q8WUW1), C3orf26 (Q9BQ75), C3orf75 (Q0PNE2), C4orf27 (Q9NWWY4), C4orf41 (Q7Z392), C5orf32 (Q9H1C7), C6orf130 (Q9Y530), C6orf211 (Q9H993), C7orf25 (Q9BPX7), C7orf28B (P86790), C7orf41 (Q8N3F0), C7orf59 (Q0VGL1), C9orf142 (Q9BUH6), C9orf23 (Q8N5L8), C9orf41 (Q8N4J0), C9orf64 (Q5T6V5), CA11 (O75493), CAB39 (Q9Y376), CACNA2D1 (P54289), CACYBP (Q9HB71), CAD (P27708), CADM1 (Q9BY67), CADM4 (Q8NFZ8), CALB1 (P05937), CALD1 (Q05682), CALM1 (P62158), CAMK2D (Q13557), CAND1 (Q86VP6), CAP1 (Q01518), CAPN1 (P07384), CAPN2 (P17655), CAPN5 (O15484), CAPNS1 (P04632), CAPS (Q13938), CAPZA1 (P52907), CAPZA2 (P47755), CAPZB (P47756),

EP 3 033 418 B1

(continued)

CARHSP1 (Q9Y2V2), CARKD (Q8IW45), CARM1 (Q86X55), CARS (P49589), CASK (O14936), CASP3 (P42574), CASP6 (P55212), CAT (P04040), CBF1 (Q13951), CBR1 (P16152), CBR3 (O75828), CBS (P35520), CBWD2 (Q8IU1), CBX1 (P83916), CBX3 (Q13185), CBX5 (P45973), CC2D1A (Q6P1 N0), CC2D1B (Q5T0F9), CCAR1 (Q8IX12), CCBL1 (Q16773), CCBL2 (Q6YP21), CCDC22 (O60826), CCDC25 (Q86WR0), CCDC53 (Q9Y3C0), CCDC56 (Q9Y2R0), CCDC93 (Q567U6), CCNC (P24863), CCND2 (P30279), CCNH (P51946), CCT2 (P78371), CCT3 (P49368), CCT4 (P50991), CCT5 (P48643), CCT6A (P40227), CCT7 (Q99832), CCT8 (P50990), CD109 (Q6YHK3), CD151 (P48509), CD276 (Q5ZPR3), CD44 (P16070), CD47 (Q08722), CD59 (P13987), CD63 (P08962), CD81 (P60033), CD9 (P21926), CD99 (P14209), CDC16 (Q13042), CDC23 (Q9UJX2), CDC27 (P30260), CDC34 (P49427), CDC37 (Q16543), CDC40 (O60508), CDC42 (P60953), CDC5L (Q99459), CDCP1 (Q9H5V8), CDH2 (P19022), CDK1 (P06493), CDK2 (P24941), CDK2AP2 (O75956), CDK4 (P11802), CDK5 (Q00535), CDK5RAP3 (Q96JB5), CDK7 (P50613), CDKN2A (P42771), CDKN2AIP (Q9NXV6), CELSR1 (Q9NYQ6), CELSR2 (Q9HCU4), CEP57 (Q86XR8), CFL1 (P23528), CFL2 (Q9Y281), CHAC2 (Q8WUX2), CHAF1B (Q13112), CHD4 (Q14839), CHEK2 (O96017), CHERP (Q8IWX8), CHID1 (Q9BWS9), CHML (P26374), CHMP1B (Q7LBR1), CHMP2A (O43633), CHMP4A (Q9BY43), CHMP4B (Q9H444), CHMP6 (Q96FZ7), CHORDC1 (Q9UHD1), CHP (Q99653), CHRAC1 (Q9NRG0), CHST14 (Q8NCH0), CHST3 (Q7LGC8), CHURC1 (Q8WUH1), CIAO1 (O76071), CIAPIN1 (Q6FI81), CIRH1A (Q969X6), CKAP5 (Q14008), CKB (P12277), CLASP1 (Q7Z460), CLDN3 (O15551), CLEC18B (Q6UXF7), CLIC1 (O00299), CLIC4 (Q9Y696), CLLD6 (Q5W111), CLNS1A (P54105), CLP1 (Q92989), CLPB (Q9H078), CLTA (P09496), CLTC (Q00610), CLU (P10909), CMAS (Q8NFW8), CMBL (Q96DG6), CMPK1 (P30085), CNBP (P62633), CNDP2 (Q96KP4), CNN2 (Q99439), CNN3 (Q15417), CNOT1 (A5YKK6), CNOT10 (Q9H9A5), CNOT6L (Q96LI5), CNOT7 (Q9UIV1), CNP (P09543), COASY (Q13057), COBRA1 (Q8WX92), COG1 (Q8WTW3), COG2 (Q14746), COG3 (Q96JB2), COG4 (Q9H9E3), COG5 (Q9UP83), COG6 (Q9Y2V7), COG7 (P83436), COG8 (Q96MW5), COL11A1 (P12107), COL14A1 (Q05707), COL6A1 (P12109), COMMD1 (Q8N668), COMMD10 (Q9Y6G5), COMMD2 (Q86X83), COMMD3 (Q9UBI1), COMMD4 (Q9H0A8), COMMD5 (Q9GZQ3), COMMD6 (Q7Z4G1), COMMD7 (Q86VX2), COMMD8 (Q9NX08), COMMD9 (Q9P000), COMT (P21964), COPA (P53621), COPB1 (P53618), COPB2 (P35606), COPE (O14579), COPG (Q9Y678), COPG2 (Q9UBF2), COPS2 (P61201), COPS3 (Q9UNS2), COPS4 (Q9BT78), COPS5 (Q92905), COPS6 (Q7L5N1), COPS7A (Q9UBW8), COPS7B (Q9H9Q2), COPS8 (Q99627), COPZ1 (P61923), CORO1A (P31146), CORO1B (Q9BR76), CORO1C (Q9ULV4), CORO2B (Q9UQ03), CORO7 (P57737), COTL1 (Q14019), COX5A (P20674), COX5B (P10606), COX6C (P09669), COX7A2 (P14406), CP (P00450), CPD (O75976), CPN2 (P22792), CPNE1 (Q99829), CPNE3 (O75131), CPNE7 (Q9UBL6), CPSF1 (Q10570), CPSF2 (Q9P2I0), CPSF3 (Q9UKF6), CPSF7 (Q8N684), CPXM1 (Q96SM3), CRIP2 (P52943), CRK (P46108),

(continued)

CRLF3 (Q8IUI8), CRTAP (O75718), CRYAB (P02511), CRYM (Q14894), CRYZ (Q08257), CRYZL1 (O95825),
 CS (O75390), CSDE1 (O75534), CSE1L (P55060), CSK (P41240), CSNK1A1 (P48729), CSNK2A1 (P68400),
 CSNK2B (P67870), CSRP1 (P21291), CSRP2 (Q16527), CSTB (P04080), CSTF1 (Q05048), CSTF2T (Q9H0L4),
 CSTF3 (Q12996), CTBP1 (Q13363), CTBP2 (P56545), CTNNA1 (P35221), CTNNB1 (P35222), CTNNBL1
 (Q8WYA6), CTNND1 (O60716), CTPS (P17812), CTPS2 (Q9NRF8), CTR9 (Q6PD62), CTSC (P53634), CTSD
 (P07339), CTSF (Q9UBX1), CTSL2 (O60911), CTU1 (Q7Z7A3), CTU2 (Q2VPK5), CUL1 (Q13616), CUL2
 (Q13617), CUL3 (Q13618), CUL4A (Q13619), CUL4B (Q13620), CUL5 (Q93034), CWF19L1 (Q69YN2), CXADR
 (P78310), CXorf26 (Q9BVG4), CYB5A (P00167), CYCS (P99999), CYFIP1 (Q7L576), CYFIP2 (Q96F07), CYR61
 (O00622), DAG1 (Q14118), DAK (Q3LXA3), DARS (P14868), DAZAP1 (Q96EP5), DBI (P07108), DBN1 (Q16643),
 DBNL (Q9UJU6), DBR1 (Q9UK59), DCAF7 (P61962), DCAF8 (Q5TAQ9), DCD (P81605), DCK (P27707), DCLK1
 (O15075), DCPS (Q96C86), DCTD (P32321), DCTN1 (Q14203), DCTN2 (Q13561), DCTN3 (O75935), DCTN4
 (Q9UJW0), DCTN5 (Q9BTE1), DCTN6 (O00399), DCUN1D1 (Q96GG9), DCUN1D5 (Q9BTE7), DCXR (Q7Z4W1),
 DDA1 (Q9BW61), DDAH2 (O95865), DDB1 (Q16531), DDB2 (Q92466), DDI2 (Q5TDH0), DDR1 (Q08345), DDT
 (P30046), DDX1 (Q92499), DDX17 (Q92841), DDX19A (Q9NUU7), DDX21 (Q9NR30), DDX23 (Q9BUQ8), DDX39
 (O00148), DDX3X (O00571), DDX5 (P17844), DDX51 (Q8N8A6), DDX6 (P26196), DECR1 (Q16698), DEF
 (Q68CQ4), DEFA1 (P59665), DENR (O43583), DERA (Q9Y315), DFFA (O00273), DHFR (P00374), DHPS
 (P49366), DHRS1 (Q96LJ7), DHRS11 (Q6UWP2), DHRS4 (Q9BTZ2), DHX15 (O43143), DHX16 (O60231),
 DHX29 (Q7Z478), DHX36 (Q9H2U1), DHX9 (Q08211), DIAPH1 (O60610), DIAPH2 (O60879), DIMIT1L
 (Q9UNQ2), DIP2B (Q9P265), DIP2C (Q9Y2E4), DIS3 (Q9Y2L1), DIS3L2 (Q8IYB7), DKC1 (O60832), DLG1
 (Q12959), DNAH17 (Q9UFH2), DNAJA1 (P31689), DNAJA2 (O60884), DNAJB1 (P25685), DNAJB4 (Q9UDY4),
 DNAJC13 (O75165), DNAJC3 (Q13217), DNAJC7 (Q99615), DNASE1L1 (P49184), DNM1 (Q05193), DNM1L
 (O00429), DNM2 (P50570), DNPEP (Q9ULA0), DOCK1 (Q14185), DOCK4 (Q8N1I0), DOCK5 (Q9H7D0), DOCK7
 (Q96N67), DOHH (Q9BU89), DOM3Z (O77932), DPCD (Q9BVM2), DPH1 (Q9BZG8), DPH2 (Q9BQC3), DPH5
 (Q9H2P9), DPM1 (O60762), DPP3 (Q9NY33), DPP9 (Q86TI2), DPY30 (Q9C005), DPYSL2 (Q16555), DPYSL3
 (Q14195), DPYSL4 (O14531), DPYSL5 (Q9BPU6), DRG1 (Q9Y295), DRG2 (P55039), DSG1 (Q02413), DSP
 (P15924), DST (Q03001), DSTN (P60981), DTD1 (Q8TEA8), DTYMK (P23919), DUS2L (Q9NX74), DUSP12
 (Q9UNI6), DUSP23 (Q9BVJ7), DUSP3 (P51452), DYM (Q7RTS9), DYNC1H1 (Q14204), DYNC1I2 (Q13409),
 DYNC1LI1 (Q9Y6G9), DYNC1LI2 (O43237), DYNC2H1 (Q8NCM8), DYNLL1 (P63167), DYNLL2 (Q96FJ2),
 DYNLRB1 (Q9NP97), DYNLT1 (P63172), ECHDC1 (Q9NTX5), ECHDC3 (Q96DC8), ECHS1 (P30084), ECM29
 (Q5VYK3), EDC4 (Q6P2E9), EEA1 (Q15075), EEF1A1 (P68104), EEF1B2 (P24534), EEF1D (P29692), EEF1E1
 (O43324), EEF1G (P26641), EEF2 (P13639), EEFSEC (P57772),

EP 3 033 418 B1

(continued)

EFEMP2 (Q95967), EFHD2 (Q96C19), EFN2 (P52799), EFTUD1 (Q7Z2Z2), EFTUD2 (Q15029), EGFR (P00533), EHD1 (Q9H4M9), EHD2 (Q9NZN4), EHD4 (Q9H223), EIF1 (P41567), EIF1AX (P47813), EIF2A (Q9BY44), EIF2AK2 (P19525), EIF2B1 (Q14232), EIF2B2 (P49770), EIF2B3 (Q9NR50), EIF2B4 (Q9UI10), EIF2B5 (Q13144), EIF2C2 (Q9UKV8), EIF2S1 (P05198), EIF2S2 (P20042), EIF2S3 (P41091), EIF3A (Q14152), EIF3B (P55884), EIF3C (Q99613), EIF3D (O15371), EIF3E (P60228), EIF3F (O00303), EIF3G (O75821), EIF3H (O15372), EIF3I (Q13347), EIF3J (O75822), EIF3K (Q9UBQ5), EIF3L (Q9Y262), EIF3M (Q7L2H7), EIF4A1 (P60842), EIF4A2 (Q14240), EIF4A3 (P38919), EIF4E (P06730), EIF4E2 (O60573), EIF4G1 (Q04637), EIF4G2 (P78344), EIF4G3 (O43432), EIF4H (Q15056), EIF5 (P55010), EIF5A (P63241), EIF5B (O60841), EIF6 (P56537), ELAC2 (Q9BQ52), ELAVL1 (Q15717), ELMO2 (Q96JJ3), ELP2 (Q6IA86), ELP3 (Q9H9T3), EMG1 (Q92979), EMILIN1 (Q9Y6C2), EML1 (O00423), EML2 (O95834), EML3 (Q32P44), EML4 (Q9HC35), ENAH (Q8N8S7), ENO1 (P06733), ENO2 (P09104), ENOPH1 (Q9UHY7), ENY2 (Q9NPA8), EPB41L2 (O43491), EPB41L3 (Q9Y2J2), EPHA2 (P29317), EPHB3 (P54753), EPHX1 (P07099), EPM2AIP1 (Q7L775), EPRS (P07814), ERH (P84090), ER11 (Q8IV48), ERI3 (O43414), ERP44 (Q9BS26), ESD (P10768), ESYT1 (Q9BSJ8), ETF1 (P62495), ETFA (P13804), ETFB (P38117), EXOC1 (Q9NV70), EXOC2 (Q96KP1), EXOC3 (O60645), EXOC4 (Q96A65), EXOC5 (O00471), EXOC6 (Q8TAG9), EXOC7 (Q9UPT5), EXOC8 (Q8IYI6), EXOSC1 (Q9Y3B2), EXOSC2 (Q13868), EXOSC3 (Q9NQT5), EXOSC4 (Q9NPD3), EXOSC5 (Q9NQT4), EXOSC6 (Q5RKV6), EXOSC7 (Q15024), EXOSC8 (Q96B26), EXOSC9 (Q06265), EXTL3 (Q43909), EYA3 (Q99504), EZR (P15311), F3 (P13726), F8 (P00451), F8A1 (P23610), FABP5 (Q01469), FABP7 (O15540), FADD (Q13158), FAF1 (Q9UNN5), FAH (P16930), FAHD2A (Q96GK7), FAM114A2 (Q9NRY5), FAM115A (Q9Y4C2), FAM120A (Q9NZB2), FAM125A (Q96EY5), FAM127A (A6ZKI3), FAM129B (Q96TA1), FAM136A (Q96C01), FAM168A (Q92567), FAM175B (Q15018), FAM188A (Q9H8M7), FAM3A (P98173), FAM3C (Q92520), FAM45B (Q6NSW5), FAM49B (Q9NUQ9), FAM82B (Q96DB5), FAM84B (Q96KN1), FAM98A (Q8NCA5), FAM98B (Q52LJ0), FARP1 (Q9Y4F1), FARP2 (O94887), FARSA (Q9Y285), FARSB (Q9NSD9), FASN (P49327), FAT1 (Q14517), FBL (P22087), FBLN2 (P98095), FBN1 (P35555), FBN2 (P35556), FBXL18 (Q96ME1), FBXO21 (O94952), FBXO22 (Q8NEZ5), FDFT1 (P37268), FDPS (P14324), FEN1 (P39748), FERMT1 (Q9BQL6), FERMT2 (Q96AC1), FGF1 (P05230), FGFRL1 (Q8N441), FGGY (Q96C11), FH (P07954), FHL1 (Q13642), FHL2 (Q14192), FHL3 (Q13643), FIS1 (Q9Y3D6), FKBP1A (P62942), FKBP3 (Q00688), FKBP4 (Q02790), FKBP5 (Q13451), FLII (Q13045), FLNA (P21333), FLNB (O75369), FLNC (Q14315), FLOT1 (O75955), FMNL2 (Q96PY5), FN3K (Q9H479), FN3KRP (Q9HA64), FNTA (P49354), FNTB (P49356), FOLR1 (P15328), FREM2 (Q5SZK8), FRMD8 (Q9BZ67), FSCN1 (Q16658), FSD1 (Q9BTV5), FTH1 (P02794), FTL (P02792), FTO (Q9C0B1), FTSJD2 (Q8N1G2), FUBP1 (Q96AE4), FUCA2 (Q9BTY2), FUK (Q8N0W3), FXR1 (P51114), G3BP1 (Q13283), G3BP2 (Q9UN86), G6PD (P11413),

EP 3 033 418 B1

(continued)

GAA (P10253), GALK1 (P51570), GALK2 (Q01415), GALNT1 (Q10472), GALNT2 (Q10471), GANAB (Q14697), GAP43 (P17677), GAPDH (P04406), GAPVD1 (Q14C86), GAR1 (Q9NY12), GARS (P41250), GART (P22102), GATSL2 (A6NHX0), GBA (P04062), GBE1 (Q04446), GCLM (P48507), GCN1L1 (Q92616), GDI1 (P31150), GDI2 (P50395), GEMIN5 (Q8TEQ6), GEMIN6 (Q8WXD5), GET4 (Q7L5D6), GFAP (P14136), GFPT1 (Q06210), GFPT2 (Q94808), GGCT (O75228), GGPS1 (Q95749), GINS1 (Q14691), GINS4 (Q9BRT9), GIPC1 (O14908), GIT1 (Q9Y2X7), GLA (P06280), GLB1 (P16278), GLB1L2 (Q8IW92), GLG1 (Q92896), GLIPR2 (Q9H4G4), GLMN (Q92990), GLO1 (Q04760), GLOD4 (Q9HC38), GLRX (P35754), GLRX3 (O76003), GLT25D1 (Q8NBJS), GLTP (Q9NZD2), GLTPD1 (Q5TA50), GLUD1 (P00367), GLUL (P15104), GMDS (O60547), GMFB (P60983), GMPPA (Q96IJ6), GMPPB (Q9Y5P6), GMPR (P36959), GMPR2 (Q9P2T1), GMPS (P49915), GNA11 (P29992), GNA13 (O14344), GNAI2 (P04899), GNAI3 (P08754), GNAQ (P50148), GNAS (Q5JWF2), GNB1 (P62873), GNB2 (P62879), GNB2L1 (P63244), GNB4 (Q9HAV0), GNE (Q9Y223), GNG12 (Q9UBI6), GNG4 (P50150), GNG5 (P63218), GNPDA1 (P46926), GNPAT1 (Q96EK6), GOLGA7 (Q7Z5G4), GOLGB1 (Q14789), GOLIM4 (O00461), GOLM1 (Q8NBJS), GOLPH3 (Q9H4A6), GORASP2 (Q9H8Y8), GPC1 (P35052), GPC4 (O75487), GPC6 (Q9Y625), GPD1L (Q8N335), GPI (P06744), GPLD1 (P80108), GPM6A (P51674), GPM6B (Q13491), GPN1 (Q9HCN4), GPR56 (Q9Y653), GPS1 (Q13098), GPX1 (P07203), GPX4 (P36969), GRB2 (P62993), GRHRP (Q9UBQ7), GRP (Q3ZCW2), GRPEL1 (Q9HAV7), GRWD1 (Q9BQ67), GSK3A (P49840), GSK3B (P49841), GSN (P06396), GSPT1 (P15170), GSS (P48637), GSTK1 (Q9Y2Q3), GSTM2 (P28161), GSTM3 (P21266), GSTM4 (Q03013), GSTO1 (P78417), GSTP1 (P09211), GSTT2 (P0CG29), GSTZ1 (O43708), GTF2F2 (P13984), GTF2H2 (Q13888), GTF2I (P78347), GTF3C1 (Q12789), GTF3C2 (Q8WUA4), GTF3C4 (Q9UKN8), GTPBP1 (O00178), GUK1 (Q16774), GYG1 (P46976), GYS1 (P13807), H2AFY (O75367), H2AFZ (P0C0S5), HADH (Q16836), HAGH (Q16775), HARS (P12081), HAT1 (O14929), HAUS3 (Q68CZ6), HAUS4 (Q9H6D7), HBA1 (P69905), HBB (P68871), HCFC1 (P51610), HDAC1 (Q13547), HDAC2 (Q92769), HDAC3 (O15379), HDHD2 (Q9H0R4), HDLBP (Q00341), HEATR1 (Q9H583), HEATR2 (Q86Y56), HEBP1 (Q9NRV9), HECTD3 (Q5T447), HEG1 (Q9ULI3), HELZ (P42694), HERC4 (Q5GLZ8), HEXB (P07686), HGS (O14964), HHIP (Q96QV1), HIBCH (Q6NVY1), HIF1AN (Q9NWT6), HINT1 (P49773), HIP1R (O75146), HIST1H1B (P16401), HIST1H1C (P16403), HIST1H2BM (Q99879), HIST1H2BO (P23527), HIST1H4A (P62805), HIST2H2AA3 (Q6FI13), HIST2H3A (Q71DI3), HK1 (P19367), HK2 (P52789), HLA-A (P30443), HLA-A (P01892), HLCS (P50747), HMGA1 (P17096), HMGB1 (P09429), HMGCL (P35914), HMGCS1 (Q01581), HMGN2 (P05204), HNRNPA1 (P09651), HNRNPA2B1 (P22626), HNRNPA3 (P51991), HNRNPAB (Q99729), HNRNPC (P07910), HNRNPD (Q14103), HNRNPF (P52597), HNRNPH1 (P31943), HNRNPH2 (P55795), HNRNPH3 (P31942), HNRNPK (P61978), HNRNPL (P14866), HNRNPM (P52272), HNRNPR (O43390),

(continued)

HNRNPU (Q00839), HNRNPUL2 (Q1KMD3), HNRPDL (O14979), HNRPLL (Q8WW9), HOOK3 (Q86VS8), HP
 (P00738), HP1BP3 (Q5SSJ5), HPCAL1 (P37235), HPRT1 (P00492), HPX (P02790), HRAS (P01112), HS6ST2
 (Q96MM7), HSD17B10 (Q99714), HSD17B4 (P51659), HSP90AA1 (P07900), HSP90AB1 (P08238), HSP90B1
 (P14625), HSPA12A (O43301), HSPA14 (Q0VDF9), HSPA1A (P08107), HSPA2 (P54652), HSPA4 (P34932),
 HSPA4L (Q95757), HSPA5 (P11021), HSPA8 (P11142), HSPA9 (P38646), HSPB1 (P04792), HSPB11 (Q9Y547),
 HSPBP1 (Q9NZL4), HSPD1 (P10809), HSPE1 (P61604), HSPG2 (P98160), HSPH1 (Q92598), HTATIP2
 (Q9BUP3), HTRA1 (Q92743), HTT (P42858), HUWE1 (Q7Z6Z7), HYOU1 (Q9Y4L1), IARS (P41252), ICAM1
 (P05362), IDE (P14735), IDH1 (O75874), IDH2 (P48735), IDI1 (Q13907), IDUA (P35475), IFI16 (Q16666), IFI35
 (P80217), IFIT5 (Q13325), IFITM3 (Q01628), IGF1R (P08069), IGF2BP2 (Q9Y6M1), IGF2BP3 (O00425), IGF2R
 (P11717), IGFBP3 (P17936), IGSF3 (O75054), IGSF8 (Q969P0), IKBKAP (O95163), IL1RAP (Q9NPH3), ILF2
 (Q12905), ILF3 (Q12906), ILK (Q13418), ILKAP (Q9H0C8), IMP4 (Q96G21), IMPA1 (P29218), IMPA2 (O14732),
 IMPAD1 (Q9NX62), IMPDH2 (P12268), INF2 (Q27J81), INPP1 (P49441), INPPL1 (O15357), INTS1 (Q8N201),
 INTS10 (Q9NVR2), INTS3 (Q68E01), INTS5 (Q6P9B9), IPO11 (Q9UI26), IPO13 (O94829), IPO4 (Q8TEX9), IPO5
 (O00410), IPO7 (O95373), IPO8 (O15397), IPO9 (Q96P70), IQGAP1 (P46940), IRF2BP2 (Q7Z5L9), IRF3
 (Q14653), IRGQ (Q8WZA9), ISG15 (P05161), ISOC1 (Q96CN7), ISPD (A4D126), ISYNA1 (Q9NPH2), ITFG3
 (Q9H0X4), ITGA2 (P17301), ITGA3 (P26006), ITGA4 (P13612), ITGA5 (P08648), ITGA6 (P23229), ITGA7
 (O13683), ITGAV (P06756), ITGB1 (P05556), ITGB4 (P16144), ITGB8 (P26012), ITPA (Q9BY32), JAM3
 (Q9BX67), JUP (P14923), KARS (Q15046), KBTBD4 (Q9NVX7), KBTBD6 (Q86V97), KCTD12 (Q96CX2), KDM1A
 (O60341), KEAP1 (O14145), KHDRBS1 (Q07666), KHSRP (Q92945), KIAA0174 (P53990), KIAA0196 (O12768),
 KIAA0319L (Q8IZA0), KIAA0664 (O75153), KIAA0776 (O94874), KIAA1033 (Q2M389), KIAA1279 (Q96EK5),
 KIAA1468 (Q9P260), KIAA1598 (A0MZ66), KIAA1797 (Q5VW36), KIAA1967 (Q8N163), KIF1A (Q12756), KIF3A
 (Q9Y496), KIF5B (P33176), KIF5C (O60282), KLC1 (Q07866), KLC2 (Q9H0B6), KLC4 (Q9NSK0), KLHDC3
 (Q9BQ90), KLHL13 (Q9P2N7), KNG1 (P01042), KNTC1 (P50748), KPNA1 (P52294), KPNA2 (P52292), KPNA3
 (O00505), KPNA4 (O00629), KPNA6 (O60684), KPNB1 (Q14974), KPRP (Q5T749), KRAS (P01116), KRIT1
 (O00522), KRT13 (P13646), KRT14 (P02533), KRT71 (Q3SY84), KTN1 (Q86UP2), L1CAM (P32004), LAGE3
 (Q14657), LAMA4 (Q16363), LAMA5 (O15230), LAMB1 (P07942), LAMC1 (P11047), LAMP1 (P11279), LAMP2
 (P13473), LANCL1 (O43813), LANCL2 (Q9NS86), LAP3 (P28838), LARP1 (Q6PKG0), LARS (Q9P2J5), LASP1
 (Q14847), LCAT (P04180), LCMT1 (Q9UIC8), LDHA (P00338), LDHB (P07195), LDLR (P01130), LEFTY2
 (O00292), LEPRE1 (Q32P28), LFNG (Q8NES3), LGALS1 (P09382), LGALS3 (P17931), LGALS3BP (Q08380),
 LHFP (Q9Y693), LIMA1 (Q9UHB6), LIMS1 (P48059), LIN7C (Q9NUP9), LIPG (Q9Y5X9), LLGL1 (Q15334),
 LMCD1 (Q9NZU5), LMNA (P02545),

(continued)

LMNB1 (P20700), LOXL4 (Q96JB6), LPL (P06858), LRBA (P50851), LRCH3 (Q96I18), LRG1 (P02750), LRP1
 (Q07954), LRRC20 (Q8TCA0), LRRC40 (Q9H9A6), LRRC47 (Q8N1G4), LRRC57 (Q8N9N7), LRSAM1
 (Q6UWE0), LRWD1 (Q9UFC0), LSM1 (O15116), LSM12 (Q3MHD2), LSM2 (Q9Y333), LSM3 (P62310), LSM4
 (Q9Y4Z0), LSM6 (P62312), LSM7 (Q9UK45), LSS (P48449), LTA4H (P09960), LTBP2 (Q14767), LTBP3
 (Q9NS15), LUM (P51884), LYPLA1 (O75608), LYPLA2 (O95372), LYPLAL1 (Q5VWZ2), M6PR (P20645), MACF1
 (Q9UPN3), MAD1L1 (Q9Y6D9), MAD2L1 (Q13257), MAEA (Q7L5Y9), MAGEE1 (Q9HCI5), MAGOHB (Q96A72),
 MALT1 (Q9UDY8), MAN1B1 (Q9UKM7), MAN2A1 (Q16706), MANBA (O00462), MAP1B (P46821), MAP1S
 (Q66K74), MAP2K1 (Q02750), MAP2K2 (P36507), MAP2K3 (P46734), MAP3K4 (Q9Y6R4), MAP4 (P27816),
 MAP4K4 (O95819), MAPK1 (P28482), MAPK12 (P53778), MAPK3 (P27361), MAPK9 (P45984), MAPKAPK2
 (P49137), MAPKSP1 (Q9UHA4), MAPRE1 (Q15691), MAPRE3 (Q9UPY8), MARCKS (P29966), MARCKSL1
 (P49006), MARK2 (Q7KZI7), MARS (P56192), MAT2A (P31153), MAT2B (Q9NZL9), MATR3 (P43243), MBD3
 (O95983), MBNL1 (Q9NR56), MCAM (P43121), MCAT (Q8IVS2), MCM2 (P49736), MCM3 (P25205), MCM4
 (P33991), MCM5 (P33992), MCM6 (Q14566), MCM7 (P33993), MCTS1 (Q9ULC4), MDH1 (P40925), MDH2
 (P40926), MDK (P21741), MDN1 (Q9NU22), ME1 (P48163), ME2 (P23368), MED1 (Q15648), MED16 (Q9Y2X0),
 MED17 (Q9NVC6), MED18 (Q9BUE0), MED20 (Q9H944), MED22 (Q15528), MED23 (Q9ULK4), MED27
 (Q6P2C8), MED30 (Q96HR3), MED31 (Q9Y3C7), MEMO1 (Q9Y316), MERIT40 (Q9NWW8), METAP1 (P53582),
 METAP2 (P50579), METT10D (Q86W50), METTL1 (Q9UBP6), METTL11A (Q9BV86), METTL13 (Q8N6R0),
 METTL2B (Q6P1Q9), METTL5 (Q9NRN9), MFAP2 (P55001), MFAP4 (P55083), MFGE8 (Q08431), MFI2
 (P08582), MGAT4B (Q9UQ53), MGAT5 (Q09328), MGEA5 (O60502), MICAL1 (Q8TDZ2), MIF (P14174), MIF4GD
 (A9UHW6), MINA (Q8IUF8), MINK1 (Q8N4C8), MIOS (Q9NXC5), MIS12 (Q9H081), MKLN1 (Q9UL63), MLTK
 (Q9NYL2), MMP14 (P50281), MMS19 (Q96T76), MOB2 (Q70IA6), MOBKL1B (Q9H8S9), MOBKL2A (Q96BX8),
 MOBKL3 (Q9Y3A3), MOCS2 (O96033), MON2 (Q7Z3U7), MORC2 (Q9Y6X9), MOV10 (Q9HCE1), MOXD1
 (Q6UVY6), MPI (P34949), MPP6 (Q9NZW5), MPRIIP (Q6WCQ1), MPST (P25325), MPZL1 (O95297), MRC2
 (Q9UBG0), MRI1 (Q9BV20), MRTO4 (Q9UKD2), MSH2 (P43246), MSN (P26038), MSTO1 (Q9BUK6), MTA1
 (Q13330), MTA2 (O94776), MTAP (Q13126), MTHFD1 (P11586), MTHFS (P49914), MTM1 (Q13496), MTMR1
 (Q13613), MTMR6 (Q9Y217), MTMR9 (Q96QG7), MTOR (P42345), MTPN (P58546), MTR (Q99707), MVD
 (P53602), MVK (Q03426), MVP (Q14764), MYADM (Q96S97), MYBBP1A (Q9BQG0), MYCBP (Q99417), MYD88
 (Q99836), MYH10 (P35580), MYH9 (P35579), MYL12B (O14950), MYL6 (P60660), MYO18A (Q92614), MYO1B
 (O43795), MYO1C (O00159), MYO1E (Q12965), MYO6 (Q9UM54), MYOF (Q9NZM1), MZT1 (Q08AG7), NAA10
 (P41227), NAA15 (Q9BXJ9), NAA16 (Q6N069), NAA20 (P61599), NAA30 (Q147X3), NAA38 (O95777), NAA50
 (Q9GZZ1), NACA (Q13765), NADSYN1 (Q6IA69), NAE1 (Q13564), NAGK (Q9UJ70),

EP 3 033 418 B1

(continued)

NAGLU (P54802), NAMPT (P43490), NANS (Q9NR45), NAP1L1 (P55209), NAP1L4 (Q99733), NAPA (P54920),
 NAPG (Q99747), NAPRT1 (Q6XQN6), NARS (O43776), NASP (P49321), NCAM1 (P13591), NCAPD2 (Q15021),
 NCAPG (Q9BPX3), NCBP1 (Q09161), NCBP2 (P52298), NCDN (Q9UBB6), NCKAP1 (Q9Y2A7), NCKIPSD
 (Q9NZQ3), NCL (P19338), NCS1 (P62166), NCSTN (Q92542), NDRG3 (Q9UGV2), NDRG4 (Q9ULP0), NDUFA2
 (O43678), NDUFA3 (O95167), NDUFA5 (Q16718), NDUFAB1 (O14561), NDUFS6 (O75380), NEDD4L (Q96PU5),
 NEFL (P07196), NEK9 (Q8TD19), NES (P48681), NF1 (P21359), NFIC (P08651), NFIX (O14938), NFKB2
 (Q00653), NHLRC2 (Q8NBF2), NHP2L1 (P55769), NID1 (P14543), NIP7 (Q9Y221), NIT1 (Q86X76), NIT2
 (Q9NQR4), NLE1 (Q9NVX2), NLGN4X (Q8N0W4), NLN (Q9BYT8), NMD3 (Q96D46), NME1 (P15531), NME2
 (P22392), NME3 (Q13232), NME7 (Q9Y5B8), NMT1 (P30419), NNMT (P40261), NOB1 (Q9ULX3), NOL11
 (Q9H8H0), NOL6 (Q9H6R4), NOMO2 (Q5JPE7), NONO (Q15233), NOP10 (Q9NPE3), NOP2 (P46087), NOTCH1
 (P46531), NOTCH3 (Q9UM47), NOVA2 (Q9UNW9), NPEPPS (P55786), NPLOC4 (Q8TAT6), NPM1 (P06748),
 NPM3 (O75607), NPTN (Q9Y639), NPW (Q8N729), NQO1 (P15559), NQO2 (P16083), NR2C2AP (Q86WQ0),
 NRAS (P01111), NRBP1 (Q9UHY1), NRBP2 (Q9NSY0), NRD1 (O43847), NRP2 (O60462), NSF (P46459),
 NSMAF (Q92636), NSMCE1 (Q8WV22), NSUN2 (Q08J23), NT5C (Q8TCD5), NT5DC1 (Q5TFE4), NTN1
 (O95631), NUBP1 (P53384), NUBP2 (Q9Y5Y2), NUCB1 (Q02818), NUDC (Q9Y266), NUDCD1 (Q96RS6),
 NUDCD2 (Q8WVJ2), NUDT1 (P36639), NUDT10 (Q8NFP7), NUDT12 (Q9BQG2), NUDT16 (Q96DE0),
 NUDT16L1 (Q9BRJ7), NUDT2 (P50583), NUDT21 (O43809), NUDT4 (Q9NZJ9), NUDT5 (Q9UKK9), NUMA1
 (Q14980), NUP188 (Q5SRE5), NUP37 (Q8NFB4), NUP43 (Q8NFB3), NUP54 (Q7Z3B4), NUP88 (Q99567),
 NUP93 (Q8N1F7), NUTF2 (P61970), NXN (Q6DKJ4), OBFC2B (Q9BQ15), OCRL (Q01968), ODZ2 (Q9NT68),
 ODZ3 (Q9P273), OGFOD1 (Q8N543), OGT (O15294), OLA1 (Q9NTK5), OLFML3 (Q9NRN5), OPA1 (O60313),
 OPLAH (O14841), OSBP (P22059), OSBPL1A (Q9BXW6), OSGEP (Q9NPF4), OTUB1 (Q96FW1), OVCA2
 (Q8WZ82), OXCT1 (P55809), OXSR1 (O95747), P4HB (P07237), PA2G4 (Q9UQ80), PAAF1 (Q9BRP4), PABPC1
 (P11940), PABPC4 (Q13310), PABPN1 (Q86U42), PACSIN2 (Q9UNF0), PACSIN3 (Q9UKS6), PAF1 (Q8N7H5),
 PAFAH1B1 (P43034), PAFAH1B2 (P68402), PAFAH1B3 (Q15102), PAICS (P22234), PAIP1 (Q9H074), PAK2
 (Q13177), PALD (Q9ULE6), PALLD (Q8WX93), PANK4 (Q9NVE7), PAPOLA (P51003), PAPSS1 (O43252), PARF
 (Q3YEC7), PARK7 (Q99497), PARN (O95453), PARP1 (P09874), PARP4 (Q9UKK3), PARVA (Q9NVD7), PBK
 (Q96KB5), PBLD (P30039), PCBP1 (Q15365), PCBP2 (Q15366), PCDHB2 (Q9Y5E7), PCDHGB4 (Q9UN71),
 PCDHGC3 (Q9UN70), PCID2 (Q5JVF3), PCMT1 (P22061), PCNA (P12004), PCOLCE2 (Q9UKZ9), PCYT2
 (Q99447), PDCD10 (Q9BUL8), PDCD2L (Q9BRP1), PDCD4 (Q53EL6), PDCD5 (O14737), PDCD6 (O75340),
 PDCD6IP (Q8WUM4), PDCL3 (Q9H2J4), PDDC1 (Q8NB37), PDE12 (Q6L8Q7), PDE6D (O43924), PDGFC
 (Q9NRA1), PDIA3 (P30101), PDIA6 (Q15084), PDLIM1 (O00151), PDLIM4 (P50479),

EP 3 033 418 B1

(continued)

PDLIM5 (Q96HC4), PDLIM7 (Q9NR12), PDRG1 (Q9NUG6), PDRO (Q6IAA8), PDS5A (Q29RF7), PDXK (O00764), PDXP (Q96GD0), PEA15 (Q15121), PEBP1 (P30086), PEF1 (Q9UBV8), PELO (Q9BRX2), PELP1 (Q8IZL8), PEPD (P12955), PFAS (O15067), PFDN2 (Q9UHV9), PFDN5 (Q99471), PFDN6 (O15212), PFKL (P17858), PFKM (P08237), PFKP (Q01813), PFN1 (P07737), PFN2 (P35080), PGAM1 (P18669), PGAM5 (Q96HS1), PGD (P52209), PGGT1B (P53609), PGK1 (P00558), PGLS (O95336), PGLYRP2 (Q96PD5), PGM1 (P36871), PGM2L1 (Q6PCE3), PGM3 (O95394), PGP (A6NDG6), PGRMC1 (O00264), PGRMC2 (O15173), PHF5A (Q7RTV0), PHGDH (O43175), PHKB (Q93100), PHLDA3 (Q9Y5J5), PHPT1 (Q9NRX4), PIK3CB (P42338), PIK3R4 (Q99570), PIN1 (Q13526), PIP4K2A (P48426), PIPOX (Q9P0Z9), PITPNB (P48739), PKM2 (P14618), PKP1 (Q13835), PLAA (Q9Y263), PLCD3 (Q8N3E9), PLCG1 (P19174), PLD3 (Q8IV08), PLEC (Q15149), PLEKHB2 (Q96CS7), PLIN3 (O60664), PLOD1 (Q02809), PLOD2 (O00469), PLOD3 (O60568), PLRG1 (O43660), PLS1 (O14651), PLS3 (P13797), PLSCR3 (Q9NRY6), PLTP (P55058), PLXNA1 (Q9UIW2), PLXNB2 (O15031), PLXND1 (Q9Y4D7), PM20D2 (Q8IYS1), PML (P29590), PMM2 (O15305), PMPCA (Q10713), PMPCB (O75439), PMVK (Q15126), PNMA2 (Q9UL42), PNO1 (Q9NRX1), PNP (P00491), PODXL (O00592), POLA1 (P09884), POLD1 (P28340), POLD2 (P49005), POLE3 (Q9NRF9), POLR1A (O95602), POLR1B (Q9H9Y6), POLR1C (O15160), POLR1D (Q9Y2S0), POLR1E (Q9GZS1), POLR2A (P24928), POLR2B (P30876), POLR2C (P19387), POLR2E (P19388), POLR2G (P62487), POLR2H (P52434), POLR2J (P52435), POLR2L (P62875), POLR3A (O14802), POLR3B (Q9NW08), POLR3C (Q9BU14), POLR3F (Q9H1D9), POP1 (Q99575), POP4 (O95707), POP5 (Q969H6), POP7 (O75817), PPA1 (Q15181), PPA2 (Q9H2U2), PPAT (Q06203), PPCS (Q9HAB8), PPIA (P62937), PPIB (P23284), PPID (Q08752), PPIF (P30405), PPIH (O43447), PPIL1 (Q9Y3C6), PPM1A (P35813), PPM1F (P49593), PPM1G (O15355), PPME1 (Q9Y570), PPP1CA (P62136), PPP1CB (P62140), PPP1CC (P36873), PPP1R7 (Q15435), PPP1R8 (Q12972), PPP2CA (P67775), PPP2CB (P62714), PPP2R1A (P30153), PPP2R2A (P63151), PPP2R4 (Q15257), PPP2R5C (O13362), PPP2R5D (O14738), PPP2R5E (O16537), PPP3CA (Q08209), PPP4C (P60510), PPP4R1 (Q8TF05), PPP5C (P53041), PPP6C (O00743), PPP6R3 (Q5H9R7), PPPDE2 (Q6ICB0), PPT1 (P50897), PPWD1 (Q96BP3), PRCP (P42785), PRDX1 (Q06830), PRDX2 (P32119), PRDX3 (P30048), PRDX5 (P30044), PRDX6 (P30041), PREP (P48147), PREPL (Q4J6C6), PRIM1 (P49642), PRIM2 (P49643), PRKACA (P17612), PRKACB (P22694), PRKAG1 (P54619), PRKAR1A (P10644), PRKAR2A (P13861), PRKAR2B (P31323), PRKDC (P78527), PRMT1 (Q99873), PRMT3 (O60678), PRMT5 (O14744), PROM1 (O43490), PROSC (O94903), PRPF19 (Q9UMS4), PRPF31 (Q8WWY3), PRPF4 (O43172), PRPF4B (O13523), PRPF8 (Q6P2Q9), PRPS1 (P60891), PRPS2 (P11908), PRPSAP1 (O14558), PRPSAP2 (O60256), PRSS23 (O95084), PRTFDC1 (Q9NRG1), PSAT1 (Q9Y617), PSMA1 (P25786), PSMA2 (P25787), PSMA3 (P25788), PSMA4 (P25789), PSMA5 (P28066), PSMA6 (P60900),

(continued)

PSMA7 (O14818), PSMB1 (P20618), PSMB2 (P49721), PSMB3 (P49720), PSMB4 (P28070), PSMB5 (P28074),
 PSMB6 (P28072), PSMB7 (Q99436), PSMB8 (P28062), PSMC1 (P62191), PSMC2 (P35998), PSMC3 (P17980),
 PSMC4 (P43686), PSMC5 (P62195), PSMC6 (P62333), PSMD1 (Q99460), PSMD10 (O75832), PSMD11
 (O00231), PSMD12 (O00232), PSMD13 (Q9UNM6), PSMD14 (O00487), PSMD2 (Q13200), PSMD3 (O43242),
 PSMD4 (P55036), PSMD5 (Q16401), PSMD6 (Q15008), PSMD7 (P51665), PSMD8 (P48556), PSMD9 (O00233),
 PSME1 (Q06323), PSME2 (Q9UL46), PSME3 (P61289), PSME4 (Q14997), PSMF1 (Q92530), PSMG1 (O95456),
 PSMG2 (Q969U7), PSMG3 (Q9BT73), PSPC1 (Q8WXF1), PSPH (P78330), PTBP1 (P26599), PTGES3 (Q15185),
 PTGFRN (Q9P2B2), PTGR1 (Q14914), PTGR2 (Q8N8N7), PTK2 (Q05397), PTK7 (Q13308), PTN (P21246),
 PTP4A1 (Q93096), PTPN1 (P18031), PTPN11 (Q06124), PTPN23 (Q9H3S7), PTPRA (P18433), PTPRG
 (P23470), PTPRZ1 (P23471), PUF60 (Q9UHX1), PUM1 (Q14671), PURB (Q96QR8), PUS7 (Q96PZ0), PVR
 (P15151), PWP1 (Q13610), PXDN (Q92626), PXK (Q7Z7A4), PYCR1 (P32322), PYCRL (Q53H96), PYGB
 (P11216), PYGL (P06737), OARS (P47897), QDPR (P09417), QKI (Q96PU8), QRIC1 (Q2TAL8), QSOX2
 (Q6ZRP7), QTRT1 (Q9BXR0), RAB10 (P61026), RAB11A (P62491), RAB11FIP1 (Q6WKZ4), RAB12 (Q6IQ22),
 RAB13 (P51153), RAB14 (P61106), RAB18 (Q9NP72), RAB1A (P62820), RAB1B (Q9H0U4), RAB21 (Q9UL25),
 RAB22A (Q9UL26), RAB23 (Q9ULC3), RAB27A (P51159), RAB2A (P61019), RAB34 (Q9BZG1), RAB35
 (Q15286), RAB3A (P20336), RAB3GAP1 (Q15042), RAB3GAP2 (Q9H2M9), RAB4A (P20338), RAB5A (P20339),
 RAB5B (P61020), RAB5C (P51148), RAB6A (P20340), RAB6B (Q9NRW1), RAB7A (P51149), RAB8A (P61006),
 RAB8B (Q92930), RABAC1 (Q9UI14), RABGAP1 (Q9Y3P9), RABGGTA (Q92696), RABGGTB (P53611), RABIF
 (P47224), RAC1 (P63000), RAD1 (O60671), RAD50 (Q92878), RAE1 (P78406), RAI14 (Q9P0K7), RALA
 (P11233), RALB (P11234), RALY (Q9UKM9), RAN (P62826), RANBP1 (P43487), RANBP2 (P49792), RANBP6
 (O60518), RANBP9 (Q96S59), RANGAP1 (P46060), RAP1A (P62834), RAP1B (P61224), RAP1GDS1 (P52306),
 RAP2B (P61225), RARS (P54136), RASA1 (P20936), RBBP4 (Q09028), RBBP5 (O15291), RBBP7 (O16576),
 RBBP9 (O75884), RBM12 (Q9NTZ6), RBM15 (Q96T37), RBM17 (Q96I25), RBM22 (Q9NW64), RBM4 (Q9BWF3),
 RBMX (P38159), RBP1 (P09455), RBPJ (Q06330), RBX1 (P62877), RCC1 (P18754), RCC2 (Q9P258), RCL
 (O43598), RDX (P35241), RECQL (P46063), REEP5 (Q00765), REEP6 (Q96HR9), REPS1 (Q96D71), RFC4
 (P35249), RFC5 (P40937), RFTN1 (Q14699), RHEB (Q15382), RHOA (P61586), RHOB (P62745), RHOC
 (P08134), RHOF (Q9HBH0), RHOG (P84095), RIC8A (Q9NPQ8), RMND5A (Q9H871), RNASEH2A (O75792),
 RNASEH2C (Q8TDP1), RNF123 (Q5XPI4), RNF20 (Q5VTR2), RNF213 (Q63HN8), RNF7 (Q9UBF6), RNGTT
 (O60942), RNH1 (P13489), RNMT (O43148), RNPEP (Q9H4A4), ROBLD3 (Q9Y2Q5), ROCK1 (Q13464), ROCK2
 (O75116), ROR1 (Q01973), RP2 (O75695), RPA1 (P27694), RPA2 (P15927), RPA3 (P35244), RPE (Q96AT9),
 RPF2 (Q9H7B2), RPL10 (P27635), RPL10A (P62906), RPL11 (P62913), RPL12 (P30050),

(continued)

RPL13 (P26373), RPL13A (P40429), RPL14 (P50914), RPL15 (P61313), RPL17 (P18621), RPL18 (Q07020),
 RPL18A (Q02543), RPL19 (P84098), RPL21 (P46778), RPL22 (P35268), RPL22L1 (Q6P5R6), RPL23 (P62829),
 RPL23A (P62750), RPL24 (P83731), RPL26 (P61254), RPL27 (P61353), RPL27A (P46776), RPL28 (P46779),
 RPL3 (P39023), RPL30 (P62888), RPL31 (P62899), RPL32 (P62910), RPL34 (P49207), RPL35 (P42766),
 RPL35A (P18077), RPL36 (Q9Y3U8), RPL36A (P83881), RPL36AL (Q969Q0), RPL37A (P61513), RPL38
 (P63173), RPL4 (P36578), RPL5 (P46777), RPL6 (Q02878), RPL7 (P18124), RPL7A (P62424), RPL8 (P62917),
 RPL9 (P32969), RPLP0 (P05388), RPLP1 (P05386), RPLP2 (P05387), RPP30 (P78346), RPP40 (Q75818),
 RPRD1A (Q96P16), RPS10 (P46783), RPS11 (P62280), RPS12 (P25398), RPS13 (P62277), RPS14 (P62263),
 RPS15 (P62841), RPS15A (P62244), RPS16 (P62249), RPS17 (P08708), RPS18 (P62269), RPS19 (P39019),
 RPS2 (P15880), RPS20 (P60866), RPS21 (P63220), RPS23 (P62266), RPS24 (P62847), RPS25 (P62851),
 RPS26 (P62854), RPS27 (P42677), RPS27A (P62979), RPS27L (Q71UM5), RPS28 (P62857), RPS29 (P62273),
 RPS3 (P23396), RPS3A (P61247), RPS4X (P62701), RPS4Y1 (P22090), RPS5 (P46782), RPS6 (P62753),
 RPS6KA3 (P51812), RPS7 (P62081), RPS8 (P62241), RPS9 (P46781), RPSA (P08865), RQCD1 (Q92600),
 Rraga (Q7L523), RRAS (P10301), RRAS2 (P62070), RRPB1 (Q9P2E9), RRM1 (P23921), RRM2 (P31350),
 RRM2B (Q7LG56), RRP12 (Q5JTH9), RRP9 (Q43818), RSL1D1 (Q76021), RSU1 (Q15404), RTCD1 (Q00442),
 RTN3 (Q95197), RTN4 (Q9NQC3), RUVBL1 (Q9Y265), RUVBL2 (Q9Y230), RWDD2B (P57060), S100A10
 (P60903), S100A11 (P31949), S100A13 (Q99584), S100A16 (Q96FQ6), S100A4 (P26447), S100A6 (P06703),
 S100A8 (P05109), SAAL1 (Q96ER3), SACS (Q9NZJ4), SAE1 (Q9UBE0), SAFB2 (Q14151), SAMHD1 (Q9Y3Z3),
 SAP18 (Q00422), SAR1A (Q9NR31), SARM1 (Q6SZW1), SARS (P49591), SART3 (Q15020), SBDS (Q9Y3A5),
 SBF1 (Q95248), SCARB1 (Q8WTV0), SCARB2 (Q14108), SCFD1 (Q8WVM8), SCLY (Q96I15), SCP2 (P22307),
 SCPEP1 (Q9HB40), SCRG1 (Q75711), SCRIB (Q14160), SCRNI (Q12765), SCRNI2 (Q96FV2), SCYL1
 (Q96KG9), SCYL2 (Q6P3W7), SDC1 (P18827), SDC2 (P34741), SDCBP (Q00560), SDF4 (Q9BRK5), SDHA
 (P31040), SDK1 (Q7Z5N4), SDSL (Q96GA7), SEC11A (P67812), SEC13 (P55735), SEC22B (Q75396), SEC23A
 (Q15436), SEC23B (Q15437), SEC231P (Q9Y6Y8), SEC24A (Q95486), SEC24B (Q95487), SEC24C (P53992),
 SEC24D (Q94855), SEC31A (Q94979), SEH1L (Q96EE3), SELH (Q8IZQ5), SEMA3A (Q14563), SEPSECS
 (Q9HD40), 40787 (Q9NVA2), 37500 (Q15019), 38596 (Q99719), 39326 (Q16181), 39692 (Q92599), 40057
 (Q9UHD8), SERBP1 (Q8NC51), SERPINA1 (P01009), SERPINA3 (P01011), SERPINA7 (P05543), SERPINB6
 (P35237), SERPINB8 (P50452), SERPINE1 (P05121), SERPINE2 (P07093), SERPING1 (P05155), SERPINH1
 (P50454), SETD3 (Q86TU7), SETD7 (Q8WTS6), SF3A1 (Q15459), SF3A2 (Q15428), SF3A3 (Q12874), SF3B1
 (Q75533), SF3B14 (Q9Y3B4), SF3B2 (Q13435), SF3B3 (Q15393), SF3B4 (Q15427), SF3B5 (Q9BWJ5), SFPQ
 (P23246), SFRP4 (Q6FHJ7), SGTA (Q43765), SH3BP4 (Q9P0V3),

(continued)

SH3GL1 (Q99961), SH3GLB1 (Q9Y371), SHBG (P04278), SHC1 (P29353), SHMT1 (P34896), SHMT2 (P34897), SHOC2 (Q9UQ13), SHPK (Q9UHHJ6), SKIV2L (Q15477), SKIV2L2 (P42285), SKP1 (P63208), SLC16A1 (P53985), SLC1A3 (P43003), SLC1A5 (Q15758), SLC29A1 (Q99808), SLC2A1 (P11166), SLC31A1 (O15431), SLC3A2 (P08195), SLC44A2 (Q8IWA5), SLC5A3 (P53794), SLC7A5 (Q01650), SLC9A3R1 (O14745), SLC9A3R2 (Q15599), SLIRP (Q9GZT3), SMAD4 (Q13485), SMARCA4 (P51532), SMARCA5 (O60264), SMARCC1 (Q92922), SMARCC2 (Q8TAQ2), SMARCD1 (Q96GM5), SMARCD2 (Q92925), SMARCE1 (Q969G3), SMC1A (Q14683), SMC2 (O95347), SMC3 (Q9UQE7), SMC4 (Q9NTJ3), SMC5 (Q8IY18), SMC6 (Q96SB8), SMCHD1 (A6NHR9), SMEK1 (Q6IN85), SMS (P52788), SMU1 (Q2TAY7), SMYD5 (Q6GMV2), SNAP23 (O00161), SNAPIN (O95295), SND1 (Q7KZF4), SNF8 (Q96H20), SNRNP200 (O75643), SNRNP40 (Q96DI7), SNRPA1 (P09661), SNRPB (P14678), SNRPD1 (P62314), SNRPD2 (P62316), SNRPD3 (P62318), SNRPE (P62304), SNRPF (P62306), SNRPG (P62308), SNTB1 (Q13884), SNUPN (O95149), SNX1 (Q13596), SNX12 (Q9UMY4), SNX17 (Q15036), SNX18 (Q96RF0), SNX2 (O60749), SNX27 (Q96L92), SNX3 (O60493), SNX5 (Q9Y5X3), SNX6 (Q9UNH7), SNX8 (Q9Y5X2), SNX9 (Q9Y5X1), SOD1 (P00441), SORD (Q00796), SORT1 (Q99523), SPAG9 (O60271), SPC24 (Q8NBT2), SPC25 (Q9HBM1), SPG21 (Q9NZD8), SPR (P35270), SPRYD4 (Q8WW59), SPTAN1 (Q13813), SPTBN1 (Q01082), SPTBN2 (O15020), SRGAP2 (O75044), SRI (P30626), SRM (P19623), SRP14 (P37108), SRP19 (P09132), SRP54 (P61011), SRP68 (Q9UHB9), SRP72 (O76094), SRP9 (P49458), SRPX (P78539), SRPX2 (O60687), SRR (Q9GZT4), SRRT (Q9BXP5), SRSF1 (Q07955), SRSF11 (Q05519), SRSF2 (Q01130), SRSF3 (P84103), SRSF6 (Q13247), SRSF7 (Q16629), SRSF9 (Q13242), SRXN1 (Q9BYN0), SSB (P05455), SSBP1 (Q04837), SSRP1 (Q08945), SSSCA1 (O60232), ST13 (P50502), STAG2 (Q8N3U4), STAM (Q92783), STAMPB (O95630), STAT1 (P42224), STAT3 (P40763), STIP1 (P31948), STK24 (Q9Y6E0), STK25 (O00506), STK38L (Q9Y2H1), STOM (P27105), STON2 (Q8WXE9), STRAP (Q9Y3F4), STUB1 (Q9UNE7), STX12 (Q86Y82), STX4 (Q12846), STX5 (Q13190), STX7 (O15400), STXBP1 (P61764), STXBP3 (O00186), STYX (Q8WUJ0), SUB1 (P53999), SUDS3 (Q9H7L9), SUGT1 (Q9Y2Z0), SUMO1 (P63165), SUPT16H (Q9Y5B9), SUPT4H1 (P63272), SUPT5H (O00267), SUPT6H (Q7KZ85), SVEP1 (Q4LDE5), SWAP70 (Q9UH65), SYMPK (Q92797), SYNCRIIP (O60506), SYNE1 (Q8NFF91), SYNE2 (Q8WXH0), SYNGR2 (O43760), SYNJ2BP (P57105), TAB1 (Q15750), TAF9 (Q9Y3D8), TAF9 (Q16594), TAGLN (Q01995), TAGLN2 (P37802), TALDO1 (P37837), TARDBP (Q13148), TARS (P26639), TATDN1 (Q6P1N9), TAX1BP3 (O14907), TBC1D13 (Q9NVG8), TBC1D15 (Q8TC07), TBC1D23 (Q9NUY8), TBC1D24 (Q9ULP9), TBC1D4 (O60343), TBC1D9B (Q66K14), TBCA (O75347), TBCB (Q99426), TBCD (Q9BTW9), TBCE (Q15813), TBL1XR1 (Q9BZK7), TCEA1 (P23193), TCEB1 (Q15369), TCEB2 (Q15370), TCERG1 (O14776), TCP1 (P17987), TDP2 (O95551), TEP1 (Q99973), TEX10 (Q9NXF1), TF (P02787), TFCEP2 (Q12800), TFG (Q92734), TFRC (P02786), TGFB1 (P01137),

(continued)

TGFB2 (P61812), TGFBI (Q15582), TGM1 (P22735), TH1L (Q8IXH7), THBS1 (P07996), THBS3 (P49746), THG1L (Q9NWX6), THOC2 (Q8NI27), THOC3 (Q96J01), THOC5 (Q13769), THOC6 (Q86W42), THOC7 (Q6I9Y2), THOP1 (P52888), THUMPD1 (Q9NXG2), THY1 (P04216), THYN1 (Q9P016), TIA1 (P31483), TIGAR (Q9NQ88), TIMM13 (Q9Y5L4), TIMM50 (Q3ZCQ8), TIMM8B (Q9Y5J9), TIMM9 (Q9Y5J7), TIMP1 (P01033), TIPRL (O75663), TKT (P29401), TLN1 (Q9Y490), TLN2 (Q9Y4G6), TM9SF2 (Q99805), TM9SF3 (Q9HD45), TMED10 (P49755), TMED2 (Q15363), TMED7 (Q9Y3B3), TMED9 (Q9BVK6), TMEM167A (Q8TBQ9), TMEM2 (Q9UHN6), TMEM50B (P56557), TMEM87A (Q8NBN3), TMOD3 (Q9NYL9), TNC (P24821), TNPO1 (Q92973), TNPO2 (O14787), TNPO3 (Q9Y5L0), TOLLIP (Q9H0E2), TOMM20 (Q15388), TOMM22 (Q9NS69), TOMM34 (Q15785), TOMM5 (Q8N4H5), TOMM70A (O94826), TOP1 (P11387), TOP2B (Q02880), TOR1B (O14657), TP53BP1 (Q12888), TP53RK (Q96S44), TPI1 (P60174), TPM3 (P06753), TPM3L (A6NL28), TPM4 (P67936), TPMT (P51580), TPP1 (O14773), TPP2 (P29144), TPR (P12270), TPRG1 L (Q5T0D9), TPRKB (Q9Y3C4), TPT1 (P13693), TRAF2 (Q12933), TRAP1 (Q12931), TRAPPC1 (Q9Y5R8), TRAPPC2L (Q9UL33), TRAPPC3 (O43617), TRAPPC4 (Q9Y296), TRAPPC5 (Q8IUR0), TRAPPC6A (O75865), TRAPPC6B (Q86SZ2), TRIM22 (Q8IYM9), TRIM25 (Q14258), TRIM28 (Q13263), TRIP12 (Q14669), TRIP13 (Q15645), TRIP6 (Q15654), TRMT1 (Q9NXH9), TRMT112 (Q9UI30), TRMT5 (Q32P41), TRMT6 (Q9UJA5), TRMT61A (Q96FX7), TRNT1 (Q96Q11), TROVE2 (P10155), TRRAP (Q9Y4A5), TSG101 (Q99816), TSKU (Q8WUA8), TSPAN14 (Q8NG11), TSPAN4 (O14817), TSPAN5 (P62079), TSPAN6 (O43657), TSPAN9 (O75954), TSSC1 (Q53HC9), TSTA3 (Q13630), TTC1 (Q99614), TTC37 (Q6PGP7), TTC38 (Q5R3I4), TTC5 (Q8N0Z6), TTC9C (Q8N5M4), TTL (Q8NG68), TTL12 (Q14166), TTN (Q8WZ42), TTR (P02766), TTYH1 (Q9H313), TTYH2 (Q9BSA4), TTYH3 (Q9C0H2), TUBA1B (P68363), TUBA1C (Q9BQE3), TUBB (P07437), TUBB2A (Q13885), TUBB2B (Q9BVA1), TUBB2C (P68371), TUBB3 (Q13509), TUBB4 (P04350), TUBB6 (Q9BUF5), TUBG1 (P23258), TUBGCP2 (Q9BSJ2), TUBGCP3 (Q96CW5), TWF1 (Q12792), TWF2 (Q6IBS0), TXN (P10599), TXNDC17 (Q9BRA2), TXNDC9 (O14530), TXNL1 (O43396), TXNL4B (Q9NX01), TXNRD1 (Q16881), TYMS (P04818), U2AF1 (Q01081), U2AF2 (P26368), UAP1 (Q16222), UBA1 (P22314), UBA2 (Q9UBT2), UBA3 (Q8TBC4), UBA5 (Q9GZZ9), UBA6 (A0AVT1), UBE2D1 (P51668), UBE2D3 (P61077), UBE2E1 (P51965), UBE2G2 (P60604), UBE2I (P63279), UBE2J2 (Q8N2K1), UBE2K (P61086), UBE2L3 (P68036), UBE2M (P61081), UBE2N (P61088), UBE2O (Q9C0C9), UBE2V1 (Q13404), UBE2V2 (Q15819), UBE2Z (Q9H832), UBE3A (Q05086), UBE4A (Q14139), UBE4B (O95155), UBL3 (O95164), UBL4A (P11441), UBL5 (Q9BZL1), UBR1 (Q8I WV7), UBR4 (Q5T4S7), UBD1 (Q9HAC8), UBXN1 (Q04323), UCHL1 (P09936), UCHL3 (P15374), UCHL5 (Q9Y5K5), UCK2 (Q9BZX2), UFC1 (Q9Y3C8), UFD1L (Q92890), UFSP2 (Q9NUQ7), UGDH (O60701), UGP2 (Q16851), UMPS (P11172), UNC119B (A6NIH7), UNC45A (Q9H3U1), UPF1 (Q92900), UPP1 (Q16831), UROD (P06132), UROS (P10746), USO1 (O60763), USP10 (Q14694), USP11 (P51784), USP14 (P54578), USP15 (Q9Y4E8), USP24 (Q9UPU5), USP39 (Q53GS9), USP5 (P45974), USP7 (Q93009), USP9X (Q93008), UTP15 (Q8TED0), UXS1 (Q8NBZ7), UXT (Q9UBK9), VAC14 (Q08AM6), VAMP3 (Q15836), VAMP5 (O95183), VAPA (Q9P0L0), VAPB (O95292), VARS (P26640), VASN (Q6EMK4), VASP (P50552), VAT1 (Q99536), VAV2 (P52735), VBP1 (P61758), VCAN (P13611), VCL (P18206), VCP (P55072), VIM (P08670), VPRBP (Q9Y4B6), VPS11 (Q9H270), VPS13C (Q709C8), VPS16 (Q9H269), VPS18 (Q9P253), VPS24 (Q9Y3E7), VPS25 (Q9BRG1), VPS26A (O75436), VPS26B (Q4G0F5), VPS28 (Q9UK41), VPS29 (Q9UBQ0), VPS33A (Q96AX1), VPS33B (Q9H267), VPS35 (Q96QK1), VPS36 (Q86VN1), VPS37B (Q9H9H4), VPS39 (Q96JC1), VPS45 (Q9NRW7), VPS4A (Q9UN37), VPS4B (O75351), VPS53 (Q5VIR6), VRK1 (Q99986), VTA1 (Q9NP79), VWA1 (Q6PCB0), VWA5A (O00534), WARS (P23381), WASF1 (Q92558), WASL (O00401), WDFY1 (Q8IWB7), WDR1 (O75083), WDR11 (Q9BZH6), WDR12 (Q9GZL7), WDR18 (Q9BV38), WDR26 (Q9H7D7), WDR33 (Q9C0J8), WDR4 (P57081), WDR43 (Q15061), WDR45L (Q5MNZ6), WDR48 (Q8TAF3), WDR5 (P61964), WDR54 (Q9H977), WDR55 (Q9H6Y2), WDR59 (Q6PJ19), WDR6 (Q9NNW5), WDR61 (Q9GZS3), WDR73 (Q6P4I2), WDR77 (Q9BQA1), WDR82 (Q6UXN9), WDR91 (A4D1P6), WDR92 (Q96MX6), WNK1 (Q9H4A3), XNPEP1 (Q9NQW7), XPO1 (O14980), XPO4 (Q9C0E2), XPO5 (Q9HAV4), XPO6 (Q96QU8), XPO7 (Q9UIA9), XPOT (O43592), XRCC1 (P18887), XRCC5 (P13010), XRCC6 (P12956), XRN2 (Q9H0D6), YARS (P54577), YBX1 (P67809), YEATS4 (O95619), YES1 (P07947), YIPF4 (Q9BSR8), YKT6 (O15498), YPEL5 (P62699), YRDC (Q86U90), YTHDF2 (Q9Y5A9), YWHAB (P31946), YWHAE (P62258), YWHAG (P61981), YWHAH (Q04917), YWHAQ (P27348), YWHAZ (P63104), ZC3HAV1L (Q96H79), ZCCHC3 (Q9NUD5), ZER1 (Q7Z7L7), ZFPL1 (O95159), ZFR (Q96KR1), ZMAT2 (Q96NC0), ZNF259 (O75312), ZW10 (O43264), ZWILCH (Q9H900), ZYG11 B (Q9C0D3), ZYX (Q15942), ZZE1 (O43149).

Table 20: 100 most abundant proteins (name and SwissProt accession number) observed in CTX0E03 exosomes

	Identified proteins	Accession number
5	Actin, cytoplasmic 2	P63261
	Glyceraldehyde-3-phosphate dehydrogenase	P04406
	Histone H4	P62805
	Pyruvate kinase isozymes M1/M2	P14618
10	Alpha-enolase	P06733
	Heat shock protein HSP 90-beta	P08238
	Ubiquitin-40S ribosomal protein S27a	P62979
	Heat shock cognate 71 kDa protein	P11142
15	Haptoglobin	P00738
	Heat shock protein HSP 90-alpha	P07900
	Phosphoglycerate kinase 1	P00558
20	Actin, alpha cardiac muscle 1	P68032
	40S ribosomal protein S3	P23396
	Elongation factor 1-alpha 1	P68104
	GTP-binding nuclear protein Ran	P62826
25	Histone H2B type 1-M	Q99879
	Peptidyl-prolyl cis-trans isomerase A	P62937
	Profilin-1	P07737
30	Elongation factor 2	P13639
	Fatty acid synthase	P49327
	Tubulin beta-2C chain	P68371
	Tubulin alpha-1 B chain	P68363
35	Tubulin beta chain	P07437
	40S ribosomal protein S11	P62280
	Eukaryotic initiation factor 4A-I	P60842
40	T-complex protein 1 subunit theta	P50990
	14-3-3 protein theta	P27348
	40S ribosomal protein S18	P62269
45	Tubulin beta-3 chain	Q13509
	T-complex protein 1 subunit beta	P78371
	40S ribosomal protein S16	P62249
	Heat shock 70 kDa protein 1A/1B	P08107
50	Histone H3.2	Q71D13
	Transketolase	P29401
	40S ribosomal protein SA	P08865
55	Clusterin	P10909
	Fatty acid-binding protein, brain	O15540
	Hemopexin	P02790

(continued)

	Identified proteins	Accession number
5	T-complex protein 1 subunit gamma	P49368
	Tubulin beta-2B chain	Q9BVA1
	Adenosylhomocysteinase	P23526
	T-complex protein 1 subunit eta	Q99832
10	40S ribosomal protein S15a	P62244
	T-complex protein 1 subunit delta	P50991
	Vimentin	P08670
15	Guanine nucleotide-binding protein subunit beta-2-like 1	P63244
	Dihydropyrimidinase-related protein 3	Q14195
	Elongation factor 1-gamma	P26641
	Fascin	Q16658
20	Creatine kinase B-type	P12277
	X-ray repair cross-complementing protein 5	P13010
	40S ribosomal protein S2	P15880
25	Histone H2A type 2-A	Q6FI13
	40S ribosomal protein S4, X isoform	P62701
	14-3-3 protein zeta/delta	P63104
	Heterogeneous nuclear ribonucleoprotein A1	P09651
30	CD81 antigen	P60033
	Keratin, type I cytoskeletal 14	P02533
	ATP-citrate synthase	P53396
35	40S ribosomal protein S9	P46781
	Transgelin-2	P37802
	Fructose-bisphosphate aldolase A	P04075
	Ubiquitin-like modifier-activating enzyme 1	P22314
40	Peroxiredoxin-1	Q06830
	40S ribosomal protein S5	P46782
	T-complex protein 1 subunit epsilon	P48643
45	60S ribosomal protein L30	P62888
	T-complex protein 1 subunit alpha	P17987
	60S ribosomal protein L12	P30050
	Cofilin-1	P23528
50	Heterogeneous nuclear ribonucleoproteins A2/B1	P22626
	Eukaryotic translation initiation factor 5A-1	P63241
	Phosphoglycerate mutase 1	P18669
55	Clathrin heavy chain 1	Q00610
	Dihydropyrimidinase-related protein 2	Q16555
	60S ribosomal protein L35a	P18077

(continued)

Identified proteins	Accession number
T-complex protein 1 subunit zeta	P40227
Carbonyl reductase [NADPH] 1	P16152
40S ribosomal protein S3a	P61247
Ferritin heavy chain	P02794
Annexin A2	P07355
Myosin light polypeptide 6	P60660
Major vault protein	Q14764
Heterogeneous nuclear ribonucleoprotein D0	Q14103
60S acidic ribosomal protein P0	P05388
X-ray repair cross-complementing protein 6	P12956
40S ribosomal protein S20	P60866
Protein arginine N-methyltransferase 1	Q99873
40S ribosomal protein S10	P46783
Transaldolase	P37837
Histone H2B type 1-	P23527
Triosephosphate isomerase	P60174
Protein S100-A6	P06703
40S ribosomal protein S17	P08708
CD9 antigen	P21926
Filamin-A	P21333
Peptidyl-prolyl cis-trans isomerase FKBP4	Q02790
Programmed cell death 6-interacting protein	Q8WUM4
Glutathione S-transferase P	P09211
14-3-3 protein epsilon	P62258

Microvesicles

[0358] 2940 proteins were identified by Mass spectrometry in Microvesicles isolated from the initial stages of an Integra culture (week 2) and purified by centrifugation at 10,000 x g. The gene names and corresponding SWISSPROT accession numbers (in brackets) of all 2940 proteins are listed in Table 21 (in alphabetical order of gene name) and the 100 most abundant proteins are listed in Table 22, in order of decreasing abundance.

Table 21: Gene names and SWISSPROT accession numbers of all 2940 proteins identified in CTX0E03 microvesicles (listed in alphabetical order of gene name).

5	A1BG (P04217), AACS (Q86V21), AAMP (Q13685), AARS (P49588), AARSD1 (Q9BTE6), AASDHPPT (Q9NRN7), ABCA3 (Q99758), ABCC1 (P33527), ABCC4 (O15439), ABCE1 (P61221), ABCF1 (Q8NE71), ABCF2 (Q9UG63), ABCF3 (Q9NUQ8), ABHD14B (Q96IU4), ABI1 (Q8IZP0), ABR (Q12979), ACAA1 (P09110), ACAA2 (P42765), ACACA (Q13085), ACADM (P11310), ACADVL (P49748), ACAT1 (P24752), ACAT2 (Q9BWD1), ACBD6 (Q9BR61), ACBD7 (Q8N6N7), ACLY (P53396), ACO1 (P21399), ACO2 (Q99798), ACOT1 (Q86TX2), ACOT13 (Q9NPJ3), ACOT7 (O00154), ACOX1 (Q15067), ACOX3 (O15254), ACP1 (P24666), ACSL1 (P33121),
10	ACSL3 (O95573), ACSL4 (O60488), ACSS2 (Q9NR19), ACTC1 (P68032), ACTG1 (P63261), ACTL6A (O96019), ACTN1 (P12814), ACTN4 (O43707), ACTR10 (Q9NZ32), ACTR1A (P61163), ACTR1B (P42025), ACTR2 (P61160), ACTR3 (P61158), ACY1 (Q03154), ADAM10 (O14672), ADAM9 (Q13443), ADAMTS15 (Q8TE58), ADAMTS16 (Q8TE57), ADAR (P55265), ADD1 (P35611), ADD3 (Q9UEY8), ADH5 (P11766), ADK (P55263), ADO (Q96SZ5), ADPRH (P54922), ADRBK1 (P25098), ADRM1 (Q16186), ADSL (P30566), ADSS (P30520),
15	AEBP1 (Q8IUX7), AFM (P43652), AGL (P35573), AGPS (O00116), AGRN (O00468), AHCY (P23526), AHCYL1 (O43865), AHNAK (Q09666), AHNAK2 (Q8IVF2), AHS1 (O95433), AHSG (P02765), AIDA (Q96BJ3), AIFM1 (O95831), AIMP1 (Q12904), AIMP2 (Q13155), AIP (O00170), AK1 (P00568), AK2 (P54819), AK3 (Q9UIJ7), AK4 (P27144), AKAP12 (Q02952), AKAP9 (Q99996), AKR1A1 (P14550), AKR1B1 (P15121), AKR1C1 (Q04828), AKR7A2 (O43488), AKR7A3 (O95154), AKT1 (P31749), ALCAM (Q13740), ALDH16A1 (Q8IZ83), ALDH18A1 (P54886),
20	ALDH2 (P05091), ALDH3A1 (P30838), ALDH7A1 (P49419), ALDH9A1 (P49189), ALDOA (P04075), ALDOC (P09972), ALKBH2 (Q6NS38), ALOX12B (O75342), AMDHD2 (Q9Y303), AMPD2 (Q01433), ANAPC1 (Q9H1A4), ANAPC4 (Q9UJX5), ANAPC5 (Q9UJX4), ANAPC7 (Q9UJX3), ANKFY1 (Q9P2R3), ANKRD17 (O75179), ANKRD28 (O15084), ANKRD52 (Q8NB46), ANP32A (P39687), ANP32B (Q92688), ANP32E (Q9BTT0),
25	ANXA1 (P04083), ANXA11 (P50995), ANXA2 (P07355), ANXA3 (P12429), ANXA4 (P09525), ANXA5 (P08758), ANXA6 (P08133), ANXA7 (P20073), AP1B1 (Q10567), AP1G1 (O43747), AP1M1 (Q9BXS5), AP1S2 (P56377), AP2A1 (O95782), AP2A2 (O94973), AP2B1 (P63010), AP2M1 (Q96CW1), AP2S1 (P53680), AP3B1 (O00203), AP3D1 (O14617), AP3M1 (Q9Y2T2), AP3S1 (Q92572), AP4S1 (Q9Y587), APEH (P13798), APEX1 (P27695), API5 (Q9BZZ5), APIP (Q96GX9), APMAP (Q9HDC9), APOA2 (P02652), APOBEC3C (Q9NRW3), APOH (P02749),
30	APOL2 (Q9BQE5), APPL1 (Q9UKG1), APRT (P07741), AQR (O60306), ARAF (P10398), ARCN1 (P48444), ARF1 (P84077), ARF4 (P18085), ARF6 (P62330), ARFGAP2 (Q8N6H7), ARFIP1 (P53367), ARFIP2 (P53365), ARG1 (P05089), ARHGAP1 (Q07960), ARHGAP5 (Q13017), ARHGDIA (P52565), ARHGEF1 (Q92888), ARHGEF10 (O15013), ARHGEF6 (Q15052), ARHGEF7 (Q14155), ARIH1 (Q9Y4X5), ARIH2 (O95376), ARL1 (P40616), ARL2 (P36404), ARL3 (P36405), ARL6IP1 (Q15041), ARL8A (Q96BM9), ARL8B (Q9NVJ2),
35	ARMC10 (Q8N2F6), ARMC6 (Q6NXE6), ARMC8 (Q8IUR7), ARMC9 (Q7Z3E5), ARPC1A (Q92747), ARPC1 B (O15143), ARPC2 (O15144), ARPC3 (O15145), ARPC4 (P59998), ARPC5 (O15511), ARPC5L (Q9BPX5), ASAH1 (Q13510), ASCC1 (Q8N9N2), ASCC3 (Q8N3C0), ASMTL (O95671), ASNA1 (O43681), ASNS (P08243), ASPSCR1 (Q9BZE9), ASS1 (P00966), ATAD3A (Q9NVI7), ATE1 (O95260), ATG101 (Q9BSB4), ATG16L1 (Q676U5), ATG3 (Q9NT62), ATG4B (Q9Y4P1), ATG7 (O95352), ATIC (P31939), ATL3 (Q6DD88), ATM (Q13315),
40	ATOX1 (O00244), ATP1A1 (P05023), ATP1B1 (P05026), ATP1B3 (P54709), ATP2A2 (P16615), ATP2B1 (P20020), ATP2B4 (P23634), ATP5A1 (P25705), ATP5B (P06576), ATP5C1 (P36542), ATP5E (P56381), ATP5F1 (P24539), ATP5H (O75947), ATP5I (P56385), ATP5L (O75964), ATP5O (P48047), ATP6AP1 (Q15904), ATP6AP2 (O75787), ATP6V0A1 (Q93050), ATP6V0D1 (P61421), ATP6V1A (P38606), ATP6V1B2 (P21281), ATP6V1C1 (P21283), ATP6V1D (Q9Y5K8), ATP6V1E1 (P36543), ATP6V1G1 (O75348), ATP6V1H (Q9UI12), ATR (Q13535),
45	ATRN (O75882), ATXN10 (Q9UBB4), B2M (P61769), B3GAT3 (O94766), B3GNT1 (O43505), BAG2 (O95816), BAG5 (Q9UL15), BAIAP2 (Q9UQB8), BANF1 (O75531), BAT1 (Q13838), BAT3 (P46379), BCAM (P50895), BCAS2 (O75934), BCAT1 (P54687), BCCIP (Q9P287), BCL2L12 (Q9HB09), BDH2 (Q9BUT1), BICD2 (Q8TD16), BLMH (Q13867), BLVRA (P53004), BLVRB (P30043), BMP1 (P13497), BOLA2 (Q9H3K6), BOP1 (Q14137), BPGM (P07738), BPNT1 (O95861), BRCC3 (P46736), BRE (Q9NXR7), BRIX1 (Q8TDN6), BROX (Q5VW32),
50	BRP16L (P0CB43), BSG (P35613), BST1 (Q10588), BTAF1 (O14981), BUB3 (O43684), BUD31 (P41223), BYSL (Q13895), BZW1 (Q7L1Q6), BZW2 (Q9Y6E2), C10orf119 (Q9BTE3), C10orf58 (Q9BRX8), C10orf76 (Q5T2E6), C11orf54 (Q9H0W9), C11orf68 (Q9H3H3), C12orf10 (Q9HB07), C12orf57 (Q99622), C14orf149 (Q96EM0), C14orf166 (Q9Y224), C14orf21 (Q86U38), C15orf58 (Q6ZNV5), C16orf13 (Q96S19), C16orf61 (Q9NRP2), C16orf80 (Q9Y6A4), C18orf21 (Q32NC0), C18orf8 (Q96DM3), C1orf123 (Q9NWW4), C1orf128 (Q9GZP4),
55	C1orf57 (Q9BSD7), C20orf11 (Q9NWU2), C20orf4 (Q9Y312),

EP 3 033 418 B1

(continued)

C21orf33 (P30042), C21orf59 (P57076), C22orf28 (Q9Y3I0), C3orf10 (Q8WUW1), C3orf26 (Q9BQ75), C3orf75 (Q0PNE2), C4orf27 (Q9NWW4), C4orf41 (Q7Z392), C4orf43 (Q96EY4), C5orf33 (Q4G0N4), C6orf211 (Q9H993), C7orf28B (P86790), C7orf50 (Q9BRJ6), C7orf59 (Q0VGL1), C8orf33 (Q9H7E9), C9orf142 (Q9BUH6), C9orf23 (Q8N5L8), C9orf41 (Q8N4J0), C9orf64 (Q5T6V5), CA11 (O75493), CA12 (O43570), CA2 (P00918), CAB39 (Q9Y376), CACNA2D1 (P54289), CACYBP (Q9HB71), CAD (P27708), CADM1 (Q9BY67), CADM4 (Q8NFZ8), CALB1 (P05937), CALD1 (Q05682), CALM1 (P62158), CALR (P27797), CALU (O43852), CAMK1 (Q14012), CAMK2D (Q13557), CAMKV (Q8NCB2), CAND1 (Q86VP6), CANX (P27824), CAP1 (Q01518), CAPN1 (P07384), CAPN2 (P17655), CAPN5 (O15484), CAPN7 (Q9Y6W3), CAPNS1 (P04632), CAPRIN1 (Q14444), CAPS (Q13938), CAPZA1 (P52907), CAPZA2 (P47755), CAPZB (P47756), CARHSP1 (Q9Y2V2), CARKD (Q8IW45), CARM1 (Q86X55), CARS (P49589), CASK (O14936), CASP14 (P31944), CASP3 (P42574), CASP7 (P55210), CAT (P04040), CBF1 (Q13951), CBR1 (P16152), CBR3 (O75828), CBS (P35520), CBX1 (P83916), CBX3 (Q13185), CBX5 (P45973), CC2D1A (Q6P1N0), CCAR1 (Q8IX12), CCBL2 (Q6YP21), CCDC102B (Q68D86), CCDC22 (O60826), CCDC25 (Q86WR0), CCDC93 (Q567U6), CCND2 (P30279), CCNY (Q8ND76), CCT2 (P78371), CCT3 (P49368), CCT4 (P50991), CCT5 (P48643), CCT6A (P40227), CCT7 (Q99832), CCT8 (P50990), CD109 (Q6YHK3), CD151 (P48509), CD276 (Q5ZPR3), CD44 (P16070), CD46 (P15529), CD47 (Q08722), CD58 (P19256), CD59 (P13987), CD63 (P08962), CD81 (P60033), CD9 (P21926), CD97 (P48960), CD99 (P14209), CDC123 (O75794), CDC16 (Q13042), CDC23 (Q9UJX2), CDC34 (P49427), CDC37 (Q16543), CDC40 (O60508), CDC42 (P60953), CDC42BPB (Q9Y5S2), CDC5L (Q99459), CDCP1 (Q9H5V8), CDH2 (P19022), CDK1 (P06493), CDK2 (P24941), CDK4 (P11802), CDK5 (Q00535), CDK5RAP3 (Q96JB5), CDK7 (P50613), CDKN2A (P42771), CDKN2AIP (Q9NXV6), CECR5 (Q9BXW7), CELF1 (Q92879), CELSR1 (Q9NYQ6), CELSR2 (Q9HCU4), CFL1 (P23528), CFL2 (Q9Y281), CHCHD3 (Q9NX63), CHD4 (Q14839), CHEK2 (O96017), CHERP (Q8IWX8), CHID1 (Q9BWS9), CHMP1A (Q9HD42), CHMP1B (Q7LBR1), CHMP2A (O43633), CHMP4A (Q9BY43), CHMP4B (Q9H444), CHMP5 (Q9NZZ3), CHMP6 (Q96FZ7), CHN1 (P15882), CHORDC1 (Q9UHD1), CHP (Q99653), CHRAC1 (Q9NRG0), CHST3 (Q7LGC8), CIAO1 (O76071), CIAPIN1 (Q6FI81), CIRBP (Q14011), CIRH1A (Q969X6), CISD2 (Q8N5K1), CKAP4 (Q07065), CKAP5 (Q14008), CKB (P12277), CLASP1 (Q7Z460), CLIC1 (O00299), CLIC4 (Q9Y696), CLLD6 (Q5W111), CLNS1A (P54105), CLPB (Q9H078), CLTA (P09496), CLTC (Q00610), CLTCL1 (P53675), CLU (P10909), CMBL (Q96DG6), CMC1 (Q7Z7K0), CMPK1 (P30085), CMTM6 (Q9NX76), CNBP (P62633), CNDP2 (Q96KP4), CNN2 (Q99439), CNN3 (Q15417), CNNM3 (Q8NE01), CNOT1 (A5YKK6), CNOT10 (Q9H9A5), CNOT6L (Q96LI5), CNP (P09543), COASY (Q13057), COBRA1 (Q8WX92), COG1 (Q8WTW3), COG3 (Q96JB2), COG4 (Q9H9E3), COG5 (Q9UP83), COG6 (Q9Y2V7), COL11A1 (P12107), COL14A1 (Q05707), COL18A1 (P39060), COL6A1 (P12109), COMMD10 (Q9Y6G5),

EP 3 033 418 B1

(continued)

COMMD2 (Q86X83), COMMD3 (Q9UBI1), COMMD5 (Q9GZQ3), COMMD8 (Q9NX08), COMMD9 (Q9P000),
 COMT (P21964), COPA (P53621), COPB1 (P53618), COPB2 (P35606), COPE (014579), COPG (Q9Y678),
 COPG2 (Q9UBF2), COPS2 (P61201), COPS3 (Q9UNS2), COPS4 (Q9BT78), COPS5 (Q92905), COPS6
 (Q7L5N1), COPS7A (Q9UBW8), COPS7B (Q9H9Q2), COPS8 (Q99627), CORO1B (Q9BR76), CORO1C
 (Q9ULV4), CORO2B (Q9UQ03), CORO7 (P57737), COTL1 (Q14019), COX4NB (O43402), COX5A (P20674),
 COX5B (P10606), COX6C (P09669), CP (P00450), CPD (O75976), CPNE1 (Q99829), CPNE2 (Q96FN4), CPNE3
 (O75131), CPNE4 (Q96A23), CPNE7 (Q9UBL6), CPOX (P36551), CPSF1 (Q10570), CPSF2 (Q9P210), CPSF3
 (Q9UKF6), CPSF3L (Q5TA45), CPSF6 (Q16630), CPSF7 (Q8N684), CPXM1 (Q96SM3), CRABP2 (P29373),
 CRIP2 (P52943), CRK (P46108), CRLF3 (Q8IUI8), CRNKL1 (Q9BZJ0), CRTAP (O75718), CRYAB (P02511),
 CRYM (Q14894), CRYZ (Q08257), CRYZL1 (O95825), CS (O75390), CSDE1 (O75534), CSE1L (P55060), CSK
 (P41240), CSNK1A1 (P48729), CSNK2A1 (P68400), CSNK2A2 (P19784), CSNK2B (P67870), CSRP1 (P21291),
 CSRP2 (Q16527), CSTB (P04080), CSTF1 (Q05048), CSTF2T (Q9H0L4), CSTF3 (Q12996), CTBP1 (Q13363),
 CTBP2 (P56545), CTNNA1 (P35221), CTNNAL1 (Q9UBT7), CTNNB1 (P35222), CTNNBL1 (Q8WYA6), CTNND1
 (O60716), CTPS (P17812), CTPS2 (Q9NRF8), CTR9 (Q6PD62), CTSC (P53634), CTSD (P07339), CTSF
 (Q9UBX1), CTSL2 (O60911), CTTN (Q14247), CTU1 (Q7Z7A3), CUL1 (Q13616), CUL2 (Q13617), CUL3
 (Q13618), CUL4A (Q13619), CUL4B (Q13620), CUL5 (Q93034), CUL7 (Q14999), CXADR (P78310), CXCL14
 (O95715), CXorf26 (Q9BVG4), CXorf38 (Q8TB03), CYB5R3 (P00387), CYC1 (P08574), CYCS (P99999), CYFIP1
 (Q7L576), CYFIP2 (Q96F07), CYR61 (O00622), DAB1 (O75553), DAD1 (P61803), DAG1 (Q14118), DAK
 (Q3LXA3), DAPK3 (O43293), DARS (P14868), DAZAP1 (Q96EP5), DBI (P07108), DBN1 (Q16643), DBNL
 (Q9UJU6), DCAF7 (P61962), DCAF8 (Q5TAQ9), DCBLD2 (Q96PD2), DCK (P27707), DCLK1 (O15075), DCPS
 (Q96C86), DCTD (P32321), DCTN1 (Q14203), DCTN2 (Q13561), DCTN3 (O75935), DCTN4 (Q9UJW0), DCTN5
 (Q9BTE1), DCTN6 (O00399), DCUN1D1 (Q96GG9), DCUN1D3 (Q8IWE4), DCUN1D5 (Q9BTE7), DCXR
 (Q7Z4W1), DDA1 (Q9BW61), DDAH1 (O94760), DDAH2 (O95865), DDB1 (Q16531), DDB2 (Q92466), DDI2
 (Q5TDH0), DDOST (P39656), DDR1 (Q08345), DDT (P30046), DDX1 (Q92499), DDX17 (Q92841), DDX18
 (Q9NVP1), DDX19A (Q9NUU7), DDX20 (Q9UHI6), DDX21 (Q9NR30), DDX23 (Q9BUQ8), DDX24 (Q9GZR7),
 DDX27 (Q96GQ7), DDX39 (O00148), DDX3X (O00571), DDX46 (Q7L014), DDX47 (Q9H0S4), DDX49 (Q9Y6V7),
 DDX5 (P17844), DDX50 (Q9BQ39), DDX51 (Q8N8A6), DDX52 (Q9Y2R4), DDX54 (Q8TDD1), DDX55 (Q8NHQ9),
 DDX56 (Q9NY93), DDX6 (P26196), DECR1 (Q16698), DECR2 (Q9NUI1), DEF (Q68CQ4), DEK (P35659), DENR
 (O43583), DERA (Q9Y315), DFFA (O00273), DFFB (O76075), DHCR24 (Q15392), DHCR7 (Q9UBM7), DHFR
 (P00374), DHPS (P49366), DHRS11 (Q6UWP2), DHRS4 (Q9BTZ2), DHX15 (O43143), DHX16 (O60231), DHX29
 (Q7Z478), DHX30 (Q7L2E3), DHX32 (Q7L7V1), DHX36 (Q9H2U1), DHX37 (Q8IY37), DHX38 (Q92620), DHX9
 (Q08211),

(continued)

DIAPH1 (O60610), DIAPH2 (O60879), DIMT1L (Q9UNQ2), DIP2A (Q14689), DIP2B (Q9P265), DIP2C (Q9Y2E4),
 DIS3 (Q9Y2L1), DIS3L2 (Q8IYB7), DKC1 (O60832), DLAT (P10515), DLD (P09622), DLG1 (Q12959), DLGAP4
 (Q9Y2H0), DLST (P36957), DMD (P11532), DNAJA1 (P31689), DNAJA2 (O60884), DNAJB1 (P25685), DNAJB11
 (Q9UBS4), DNAJB4 (Q9UDY4), DNAJB6 (O75190), DNAJC13 (O75165), DNAJC2 (Q99543), DNAJC3 (Q13217),
 DNAJC7 (Q99615), DNASE1L1 (P49184), DNM1 (Q05193), DNM1L (O00429), DNM2 (P50570), DNMT1
 (P26358), DNPEP (Q9ULA0), DOCK1 (Q14185), DOCK4 (Q8N1I0), DOCK5 (Q9H7D0), DOCK7 (Q96N67),
 DOCK9 (Q9BZ29), DOHH (Q9BU89), DPCD (Q9BVM2), DPH2 (Q9BQC3), DPH5 (Q9H2P9), DPM1 (O60762),
 DPM3 (Q9P2X0), DPP3 (Q9NY33), DPP9 (Q86TI2), DPY30 (Q9C005), DPYSL2 (Q16555), DPYSL3 (Q14195),
 DPYSL4 (O14531), DPYSL5 (Q9BPU6), DRG1 (Q9Y295), DRG2 (P55039), DSC1 (Q08554), DSG1 (Q02413),
 DSP (P15924), DST (Q03001), DSTN (P60981), DTD1 (Q8TEA8), DTNA (Q9Y4J8), DTYMK (P23919), DUS2L
 (Q9NX74), DUS3L (Q96G46), DUSP12 (Q9UNI6), DUSP3 (P51452), DYM (Q7RTS9), DYNC1H1 (Q14204),
 DYNC1I2 (Q13409), DYNC1LI1 (Q9Y6G9), DYNC1LI2 (O43237), DYNC2H1 (Q8NCM8), DYNLL1 (P63167),
 DYNLL2 (Q96FJ2), DYNLRB1 (Q9NP97), DYNLT1 (P63172), EBNA1BP2 (Q99848), ECE1 (P42892), ECHDC1
 (Q9NTX5), ECHS1 (P30084), ECM29 (Q5VYK3), EDC3 (Q96F86), EDC4 (Q6P2E9), EEA1 (Q15075), EEF1A1
 (P68104), EEF1B2 (P24534), EEF1D (P29692), EEF1E1 (O43324), EEF1G (P26641), EEF2 (P13639), EEF2K
 (O00418), EEFSEC (P57772), EFEMP2 (O95967), EFHD2 (Q96C19), EFTUD1 (Q7Z2Z2), EFTUD2 (Q15029),
 EGFR (P00533), EHD1 (Q9H4M9), EHD2 (Q9NZN4), EHD3 (Q9NZN3), EHD4 (Q9H223), EIF1AX (P47813),
 EIF2A (Q9BY44), EIF2AK2 (P19525), EIF2AK4 (Q9P2K8), EIF2B1 (Q14232), EIF2B2 (P49770), EIF2B3
 (Q9NR50), EIF2B4 (Q9UI10), EIF2B5 (Q13144), EIF2C1 (Q9UL18), EIF2C2 (Q9UKV8), EIF2S1 (P05198), EIF2S2
 (P20042), EIF2S3 (P41091), EIF3A (Q14152), EIF3B (P55884), EIF3C (Q99613), EIF3D (O15371), EIF3E
 (P60228), EIF3F (O00303), EIF3G (O75821), EIF3H (O15372), EIF3I (Q13347), EIF3J (O75822), EIF3K (Q9UBQ5),
 EIF3L (Q9Y262), EIF3M (Q7L2H7), EIF4A1 (P60842), EIF4A2 (Q14240), EIF4A3 (P38919), EIF4E (P06730),
 EIF4G1 (Q04637), EIF4G2 (P78344), EIF4H (Q15056), EIF5 (P55010), EIF5A (P63241), EIF5B (O60841), EIF6
 (P56537), ELAC2 (Q9BQ52), ELAVL1 (Q15717), ELMO2 (Q96JJ3), ELP2 (Q6IA86), ELP3 (Q9H9T3), EMD
 (P50402), EMG1 (Q92979), EML1 (O00423), EML2 (O95834), EML3 (Q32P44), EML4 (Q9HC35), ENAH
 (Q8N8S7), ENC1 (O14682), ENO1 (P06733), ENO2 (P09104), ENOPH1 (Q9UHY7), ENY2 (Q9NPA8), EPB41L2
 (O43491), EPB41L3 (Q9Y2J2), EPDR1 (Q9UM22), EPHA2 (P29317), EPHB2 (P29323), EPHB3 (P54753), EPHB4
 (P54760), EPHX1 (P07099), EPM2AIP1 (Q7L775), EPN1 (Q9Y6I3), EPRS (P07814), ERBB2IP (Q96RT1),
 ERGIC1 (Q969X5), ERH (P84090), ERI1 (Q8IV48), ERI3 (O43414), ERLIN2 (O94905), ERO1L (Q96HE7), ERP29
 (P30040), ERP44 (Q9BS26), ESD (P10768), ESYT1 (Q9BSJ8), ETF1 (P62495), ETFA (P13804), ETFB (P38117),
 EXOC1 (Q9NV70), EXOC2 (Q96KP1), EXOC3 (O60645),

EP 3 033 418 B1

(continued)

EXOC4 (Q96A65), EXOC5 (O00471), EXOC6 (Q8TAG9), EXOC6B (Q9Y2D4), EXOC7 (Q9UPT5), EXOC8 (Q8IYI6), EXOSC1 (Q9Y3B2), EXOSC10 (Q01780), EXOSC2 (Q13868), EXOSC3 (Q9NQ5T), EXOSC4 (Q9NPD3), EXOSC5 (Q9NQ5T), EXOSC6 (Q5RKV6), EXOSC7 (Q15024), EXOSC8 (Q96B26), EXOSC9 (Q06265), EZR (P15311), F11R (Q9Y624), F8 (P00451), F8A1 (P23610), FABP5 (Q01469), FABP7 (O15540), FADD (Q13158), FAH (P16930), FAHD1 (Q6P587), FAHD2A (Q96GK7), FAM115A (Q9Y4C2), FAM120A (Q9NZB2), FAM125A (Q96EY5), FAM127A (A6ZKI3), FAM129A (Q9BZQ8), FAM129B (Q96TA1), FAM136A (Q96C01), FAM175B (Q15018), FAM3C (Q92520), FAM45B (Q6NSW5), FAM49B (Q9NUQ9), FAM82B (Q96DB5), FAM84B (Q96KN1), FAM96B (Q9Y3D0), FAM98A (Q8NCA5), FAM98B (Q52LJ0), FANCI (Q9NVI1), FAR1 (Q8WVX9), FARP1 (Q9Y4F1), FARP2 (O94887), FARSA (Q9Y285), FARSB (Q9NSD9), FAS (P25445), FASN (P49327), FAT1 (Q14517), FAU (P62861), FBL (P22087), FBLN2 (P98095), FBN1 (P35555), FBN2 (P35556), FBXL18 (Q96ME1), FBXO21 (O94952), FBXO22 (Q8NEZ5), FBXW11 (Q9UKB1), FCF1 (Q9Y324), FDFT1 (P37268), FDPS (P14324), FDXR (P22570), FEN1 (P39748), FERMT1 (Q9BQL6), FERMT2 (Q96AC1), FFR (Q9UID3), FGFBP3 (Q8TAT2), FH (P07954), FHL1 (Q13642), FHL2 (Q14192), FHL3 (Q13643), FIBP (Q04327), FKBP10 (Q96AY3), FKBP1A (P62942), FKBP2 (P26885), FKBP3 (Q00688), FKBP4 (Q02790), FKBP5 (Q13451), FLG (P20930), FLG2 (Q5D862), FLII (Q13045), FLNA (P21333), FLNB (O75369), FLNC (Q14315), FLOT1 (O75955), FLOT2 (Q14254), FMNL2 (Q96PY5), FN3K (Q9H479), FN3KRP (Q9HA64), FNTA (P49354), FNTB (P49356), FOLR1 (P15328), FREM2 (Q5SZK8), FRG1 (Q14331), FRMD5 (Q7Z6J6), FRMD8 (Q9BZ67), FRYL (O94915), FSCN1 (Q16658), FSD1 (Q9BT5V), FTH1 (P02794), FTL (P02792), FTO (Q9C0B1), FTSJD2 (Q8N1G2), FUBP1 (Q96AE4), FUBP3 (Q96I24), FUCA2 (Q9BTY2), FUK (Q8N0W3), FUS (P35637), FXR1 (P51114), FXR2 (P51116), FYCO1 (Q9BQS8), FYN (P06241), G3BP1 (Q13283), G3BP2 (Q9UN86), G6PD (P11413), GAA (P10253), GALK1 (P51570), GALK2 (Q01415), GALNT1 (Q10472), GALNT2 (Q10471), GALNT7 (Q86SF2), GAN (Q9H2C0), GANAB (Q14697), GAP43 (P17677), GAPDH (P04406), GAPVD1 (Q14C86), GAR1 (Q9NY12), GARS (P41250), GART (P22102), GATSL2 (A6NHX0), GBA (P04062), GBE1 (Q04446), GBF1 (Q92538), GCDH (Q92947), GCLC (P48506), GCLM (P48507), GCN1L1 (Q92616), GDI1 (P31150), GDI2 (P50395), GEMIN4 (P57678), GEMIN5 (Q8TEQ6), GEMIN6 (Q8WXD5), GET4 (Q7L5D6), GFAP (P14136), GFM1 (Q96RP9), GFPT1 (Q06210), GFPT2 (O94808), GGCT (O75223), GGPS1 (O95749), GINS1 (Q14691), GINS2 (Q9Y248), GINS4 (Q9BRT9), GIPC1 (O14908), GIT1 (Q9Y2X7), GLA (P06280), GLB1L2 (Q8IW92), GLE1 (Q53GS7), GLG1 (Q92896), GLIPR2 (Q9H4G4), GLMN (Q92990), GLO1 (Q04760), GLOD4 (Q9HC38), GLRX (P35754), GLRX3 (O76003), GLT25D1 (Q8NBJ5), GLT25D2 (Q8IYK4), GLTP (Q9NZD2), GLUD1 (P00367), GLUL (P15104), GMDS (O60547), GMFB (P60983), GMPPA (Q96IJ6), GMPPB (Q9Y5P6), GMPR (P36959), GMPR2 (Q9P2T1), GMPS (P49915), GNA11 (P29992), GNA12 (Q03113), GNA13 (Q14344), GNAI1 (P63096),

(continued)

GNAI2 (P04899), GNAI3 (P08754), GNAQ (P50148), GNAS (Q5JWF2), GNB1 (P62873), GNB1L (Q9BYB4),
 GNB2 (P62879), GNB2L1 (P63244), GNB4 (Q9HAV0), GNE (Q9Y223), GNG10 (P50151), GNG12 (Q9UBI6),
 GNG4 (P50150), GNG5 (P63218), GNL3 (Q9BVP2), GNPDA1 (P46926), GNPAT1 (Q96EK6), GOLGA7
 (Q7Z5G4), GOLM1 (Q8NB4), GOLPH3 (Q9H4A6), GORASP2 (Q9H8Y8), GOT1 (P17174), GOT2 (P00505),
 GPC1 (P35052), GPC4 (O75487), GPC6 (Q9Y625), GPD1L (Q8N335), GPHN (Q9NQX3), GPI (P06744), GPM6A
 (P51674), GPN1 (Q9HCN4), GPR50 (Q13585), GPR56 (Q9Y653), GPS1 (Q13098), GPSM1 (Q86YR5), GPX1
 (P07203), GPX4 (P36969), GRB2 (P62993), GRHR (Q9UBQ7), GRP (Q3ZCW2), GRWD1 (Q9BQ67), GSDMA
 (Q96QA5), GSK3A (P49840), GSK3B (P49841), GSN (P06396), GSPT1 (P15170), GSR (P00390), GSS (P48637),
 GSTK1 (Q9Y2Q3), GSTM2 (P28161), GSTM3 (P21266), GSTM4 (Q03013), GSTO1 (P78417), GSTP1 (P09211),
 GSTT2 (P0CG29), GSTZ1 (O43708), GTF2E2 (P29084), GTF2F2 (P13984), GTF2H3 (Q13889), GTF21
 (P78347), GTF3C2 (Q8WUA4), GTF3C3 (Q9Y5Q9), GTF3C4 (Q9UKN8), GTPBP1 (O00178), GTPBP4 (Q9BZE4),
 GUK1 (Q16774), GYG1 (P46976), GYS1 (P13807), H1F0 (P07305), H1FX (Q92522), H2AFX (P16104), H2AFY
 (O75367), H2AFZ (P0C0S5), HADH (Q16836), HADHA (P40939), HARS (P12081), HAT1 (O14929), HAUS3
 (Q68CZ6), HAUS4 (Q9H6D7), HBA1 (P69905), HBB (P68871), HBS1L (Q9Y450), HBXIP (O43504), HCFC1
 (P51610), HDAC1 (Q13547), HDAC2 (Q92769), HDCC2 (Q7Z4H3), HDGF (P51858), HDGFRP2 (Q7Z4V5),
 HDHD2 (Q9H0R4), HDLBP (Q00341), HEATR1 (Q9H583), HEATR2 (Q86Y56), HEBP1 (Q9NRV9), HECTD3
 (Q5T447), HERC4 (Q5GLZ8), HEXB (P07686), HGS (O14964), HHIP (Q96QV1), HINT1 (P49773), HINT2
 (Q9BX68), HINT3 (Q9NQE9), HIP1R (O75146), HIST1H1B (P16401), HIST1H1C (P16403), HIST1H1D (P16402),
 HIST1H1E (P10412), HIST1H2AD (P20671), HIST1H2BJ (P06899), HIST1H2BM (Q99879), HIST1H2BO
 (P23527), HIST1H4A (P62805), HIST2H2AA3 (Q6FI13), HIST2H2AB (Q8IUE6), HIST2H2BE (Q16778),
 HIST2H3A (Q71DI3), HIST3H2BB (Q8N257), HK1 (P19367), HK2 (P52789), HLA-A (P30443), HLA-A (P01892),
 HLA-B (P03989), HMGA1 (P17096), HMGB1 (P09429), HMGB2 (P26583), HMGCL (P35914), HMGCS1
 (Q01581), HMGN1 (P05114), HMGN2 (P05204), HMGN4 (O00479), HNRNPA0 (Q13151), HNRNPA1 (P09651),
 HNRNPA2B1 (P22626), HNRNPA3 (P51991), HNRNPAB (Q99729), HNRNPC (P07910), HNRNPD (Q14103),
 HNRNPF (P52597), HNRNPH1 (P31943), HNRNPH2 (P55795), HNRNPH3 (P31942), HNRNPK (P61978),
 HNRNPL (P14866), HNRNPM (P52272), HNRNPR (O43390), HNRNPU (Q00839), HNRNPUL1 (Q9BUJ2),
 HNRNPUL2 (Q1KMD3), HNRPDL (O14979), HNRPLL (Q8WW9), HOOK3 (Q86VS8), HP (P00738), HP1BP3
 (Q5SSJ5), HPCAL1 (P37235), HPRT1 (P00492), HPX (P02790), HRAS (P01112), HRNR (Q86YZ3), HSD17B10
 (Q99714), HSD17B12 (Q53GQ0), HSD17B4 (P51659), HSDL2 (Q6YN16), HSP90AA1 (P07900), HSP90AB1
 (P08238), HSP90B1 (P14625), HSPA12A (O43301), HSPA14 (Q0VDF9), HSPA1A (P08107), HSPA4 (P34932),
 HSPA4L (Q95757), HSPA5 (P11021), HSPA8 (P11142), HSPA9 (P38646), HSPB1 (P04792), HSPBP1 (Q9NZL4),

(continued)

HSPD1 (P10809), HSPE1 (P61604), HSPG2 (P98160), HSPH1 (Q92598), HTRA1 (Q92743), HTT (P42858),
 HUWE1 (Q7Z6Z7), HYOU1 (Q9Y4L1), IARS (P41252), ICAM1 (P05362), IDE (P14735), IDH1 (O75874), IDH2
 (P48735), IDH3A (P50213), IDI1 (Q13907), IFI16 (Q16666), IFIT5 (Q13325), IFITM3 (Q01628), IFRD2 (Q12894),
 IFT172 (Q9UG01), IGF1R (P08069), IGF2BP2 (Q9Y6M1), IGF2BP3 (O00425), IGF2R (P11717), IGFBP3
 (P17936), IGFBP5 (P24593), IGHG1 (P01857), IGHG2 (P01859), IGSF3 (O75054), IGSF8 (Q969P0), IKBKAP
 (Q95163), IKBKB (O14920), IL1RAP (Q9NPH3), ILF2 (Q12905), ILF3 (Q12906), ILK (Q13418), ILKAP (Q9H0C8),
 IMMT (Q16891), IMP3 (Q9NV31), IMPA1 (P29218), IMPA2 (O14732), IMPAD1 (Q9NX62), IMPDH1 (P20839),
 IMPDH2 (P12268), INA (Q16352), INF2 (Q27J81), INPP1 (P49441), INPPL1 (O15357), INTS10 (Q9NVR2), INTS3
 (Q68E01), INTS7 (Q9NVH2), INTS8 (Q75QN2), IPO11 (Q9UI26), IPO4 (Q8TEX9), IPO5 (O00410), IPO7
 (Q95373), IPO8 (O15397), IPO9 (Q96P70), IQGAP1 (P46940), IRF2BP2 (Q7Z5L9), IRF3 (Q14653), IRGQ
 (Q8WZA9), ISOC1 (Q96CN7), ISYNA1 (Q9NPH2), ITFG3 (Q9H0X4), ITGA2 (P17301), ITGA3 (P26006), ITGA4
 (P13612), ITGA5 (P08648), ITGA6 (P23229), ITGA7 (Q13683), ITGAV (P06756), ITGB1 (P05556), ITGB1BP1
 (O14713), ITGB3 (P05106), ITGB4 (P16144), ITGB5 (P18084), ITGB8 (P26012), ITPA (Q9BY32), JAM3
 (Q9BX67), JUP (P14923), KARS (Q15046), KATNB1 (Q9BVA0), KBTBD6 (Q86V97), KCTD21 (Q4G0X4), KDM1A
 (O60341), KEAP1 (Q14145), KHDRBS1 (Q07666), KHSRP (Q92945), KIAA0020 (Q15397), KIAA0090 (Q8N766),
 KIAA0174 (P53990), KIAA0196 (Q12768), KIAA0664 (O75153), KIAA0776 (O94874), KIAA1033 (Q2M389),
 KIAA1279 (Q96EK5), KIAA1598 (A0MZ66), KIAA1797 (Q5VW36), KIAA1949 (Q6NYC8), KIAA1967 (Q8N163),
 KIDINS220 (Q9ULH0), KIF1A (Q12756), KIF2A (O00139), KIF5B (P33176), KIF5C (O60282), KLC1 (Q07866),
 KLHDC4 (Q8TBB5), KLHL13 (Q9P2N7), KLHL22 (Q53GT1), KLHL26 (Q53HC5), KNTC1 (P50748), KPNA1
 (P52294), KPNA2 (P52292), KPNA3 (O00505), KPNA4 (O00629), KPNA6 (O60684), KPNB1 (Q14974), KPRP
 (Q5T749), KRAS (P01116), KRIT1 (O00522), KRT13 (P13646), KRT14 (P02533), KRT71 (Q3SY84), KTN1
 (Q86UP2), L1CAM (P32004), LACTB2 (Q53H82), LAMA1 (P25391), LAMA4 (Q16363), LAMA5 (O15230), LAMB1
 (P07942), LAMB2 (P55268), LAMC1 (P11047), LAMP1 (P11279), LAMP2 (P13473), LANCL1 (O43813), LANCL2
 (Q9NS86), LAP3 (P28838), LARP1 (Q6PKG0), LARS (Q9P2J5), LAS1L (Q9Y4W2), LASP1 (Q14847), LBR
 (Q14739), LCMT1 (Q9UIC8), LDHA (P00338), LDHB (P07195), LDLR (P01130), LEFTY2 (O00292), LEPRE1
 (Q32P28), LGALS1 (P09382), LGALS3 (P17931), LGALS3BP (Q08380), LGALS7 (P47929), LIMA1 (Q9UHB6),
 LIMS1 (P48059), LIN7C (Q9NUP9), LIPG (Q9Y5X9), LLGL1 (Q15334), LMAN1 (P49257), LMAN2 (Q12907),
 LMCD1 (Q9NZU5), LMNA (P02545), LMNB1 (P20700), LMNB2 (Q03252), LNPEP (Q9UIQ6), LOH12CR1
 (Q969J3), LONP1 (P36776), LOR (P23490), LOXL4 (Q96JB6), LPHN2 (O95490), LPL (P06858), LRBA (P50851),
 LRG1 (P02750), LRP1 (Q07954), LRPPRC (P42704), LRRC1 (Q9BTT6), LRRC40 (Q9H9A6), LRRC47 (Q8N1G4),
 LRRC57 (Q8N9N7), LRRC59 (Q96AG4), LRRC8A (Q8IWT6), LRSAM1 (Q6UWE0),

EP 3 033 418 B1

(continued)

LSM1 (015116), LSM12 (Q3MHD2), LSM2 (Q9Y333), LSM4 (Q9Y4Z0), LSM6 (P62312), LSM7 (Q9UK45), LSS (P48449), LTA4H (P09960), LTBP2 (Q14767), LTBP3 (Q9NS15), LTN1 (O94822), LUC7L (Q9NQ29), LUC7L2 (Q9Y383), LUC7L3 (O95232), LYAR (Q9NX58), LYPLA1 (O75608), LYPLA2 (O95372), LYPLAL1 (Q5VWZ2), LZTR1 (Q8N653), M6PR (P20645), MACF1 (Q9UPN3), MACF1 (Q96PK2), MACROD1 (Q9BQ69), MAD1L1 (Q9Y6D9), MAD2L1 (Q13257), MAGEE1 (Q9HCI5), MAK16 (Q9BXY0), MALT1 (Q9UDY8), MAN1A2 (O60476), MAN1B1 (Q9UKM7), MAN2C1 (Q9NTJ4), MAP1B (P46821), MAP1LC3A (Q9H492), MAP1LC3B2 (A6NCE7), MAP2K1 (Q02750), MAP2K2 (P36507), MAP2K3 (P46734), MAP2K4 (P45985), MAP2K7 (O14733), MAP4 (P27816), MAP4K4 (O95819), MAPK1 (P28482), MAPK14 (O16539), MAPK3 (P27361), MAPKSP1 (Q9UHA4), MAPRE1 (Q15691), MAPRE3 (Q9UPY8), MARCKS (P29966), MARCKSL1 (P49006), MARK2 (Q7KZI7), MARS (P56192), MAT2A (P31153), MAT2B (Q9NZL9), MATR3 (P43243), MBD3 (Q95983), MBLAC2 (Q68D91), MBNL1 (Q9NR56), MBNL2 (Q5VZF2), MCAM (P43121), MCM2 (P49736), MCM3 (P25205), MCM4 (P33991), MCM5 (P33992), MCM6 (Q14566), MCM7 (P33993), MCTS1 (Q9ULC4), MDH1 (P40925), MDH2 (P40926), MDK (P21741), MDN1 (Q9NU22), ME1 (P48163), ME2 (P23368), MED1 (Q15648), MED10 (Q9BTT4), MED11 (Q9P086), MED17 (Q9NVC6), MED18 (Q9BUE0), MED20 (Q9H944), MED23 (Q9ULK4), MED24 (O75448), MED28 (Q9H204), MED31 (Q9Y3C7), MEMO1 (Q9Y316), MEN1 (O00255), MERIT40 (Q9NWW8), METAP1 (P53582), METAP2 (P50579), METHOD (Q86W50), METTL1 (Q9UBP6), METTL11A (Q9BV86), METTL13 (Q8N6R0), METTL2B (Q6P1Q9), METTL5 (Q9NRN9), METTL9 (Q9H1A3), MFAP2 (P55001), MFAP4 (P55083), MFGE8 (Q08431), MFI2 (P08582), MGEA5 (O60502), MICA (Q29983), MICAL1 (Q8TDZ2), MIF (P14174), MINA (Q8IUF8), MIOS (Q9NXC5), MKI67IP (Q9BYG3), MLEC (Q14165), MLLT4 (P55196), MLST8 (Q9BVC4), MLTK (Q9NYL2), MMP14 (P50281), MMP2 (P08253), MMS19 (Q96T76), MOB2 (Q70IA6), MOBKL1B (Q9H8S9), MOBKL2A (Q96BX8), MOBKL3 (Q9Y3A3), MOCS2 (O96033), MOGS (Q13724), MON2 (Q7Z3U7), MORC2 (Q9Y6X9), MOV10 (Q9HCE1), MOXD1 (Q6UVY6), MPG (P29372), MPI (P34949), MPP6 (Q9NZW5), MPRIP (Q6WCQ1), MPST (P25325), MPZL1 (O95297), MRC2 (Q9UBG0), MRE11A (P49959), MRI1 (Q9BV20), MRPS27 (Q92552), MRPS28 (Q9Y2Q9), MRPS33 (Q9Y291), MRPS34 (P82930), MRPS6 (P82932), MRTO4 (Q9UKD2), MSH2 (P43246), MSH3 (P20585), MSH6 (P52701), MSN (P26038), MSTO1 (Q9BUK6), MTA1 (Q13330), MTA2 (O94776), MTAP (Q13126), MTHFD1 (P11586), MTHFS (P49914), MTM1 (Q13496), MTMR1 (Q13613), MTMR2 (Q13614), MTMR6 (Q9Y217), MTMR9 (Q96QG7), MTOR (P42345), MTPN (P58546), MTR (Q99707), MTRR (Q9UBK8), MVD (P53602), MVK (Q03426), MVP (Q14764), MX1 (P20591), MYADM (Q96S97), MYBBP1A (Q9BQG0), MYCBP (Q99417), MYD88 (Q99836), MYH10 (P35580), MYH14 (Q7Z406), MYH9 (P35579), MYL12B (O14950), MYL6 (P60660), MYO18A (Q92614), MYO1B (O43795), MYO1C (O00159), MYO1E (Q12965), MYO5A (Q9Y4I1), MYO6 (Q9UM54), MYOF (Q9NZM1), NAA10 (P41227), NAA15 (Q9BXJ9),

(continued)

NAA16 (Q6N069), NAA25 (Q14CX7), NAA38 (Q95777), NAA50 (Q9GZZ1), NACA (Q13765), NAE1 (Q13564),
 NAGK (Q9UJ70), NAGLU (P54802), NAMPT (P43490), NANS (Q9NR45), NAP1L1 (P55209), NAP1L4 (Q99733),
 NAPA (P54920), NAPG (Q99747), NAPRT1 (Q6XQN6), NARFL (Q9H6Q4), NARS (O43776), NASP (P49321),
 NAT10 (Q9H0A0), NAT9 (Q9BTE0), NCAM1 (P13591), NCAN (O14594), NCAPD2 (Q15021), NCAPG (Q9BPX3),
 NCBP1 (Q09161), NCCRP1 (Q6ZVX7), NCDN (Q9UBB6), NCKAP1 (Q9Y2A7), NCKIPSD (Q9NZQ3), NCL
 (P19338), NCLN (Q969V3), NCS1 (P62166), NCSTN (Q92542), NDOR1 (Q9UHB4), NDRG3 (Q9UGV2), NDRG4
 (Q9ULP0), NDUFA2 (O43678), NDUFA7 (O95182), NDUFAB1 (O14561), NDUFB4 (O95168), NDUFC2 (O95298),
 NDUFS5 (O43920), NDUFS6 (O75380), NEDD8 (Q15843), NEFL (P07196), NEFM (P07197), NEK6 (Q9HC98),
 NEK9 (Q8TD19), NES (P48681), NF1 (P21359), NF2 (P35240), NFIX (Q14938), NHLRC2 (Q8NBF2), NHP2L1
 (P55769), NID1 (P14543), NIP7 (Q9Y221), NIPSNAP1 (Q9BPW8), NIT1 (Q86X76), NIT2 (Q9NQR4), NKRF
 (O15226), NLE1 (Q9NVX2), NLGN4X (Q8N0W4), NLN (Q9BYT8), NMD3 (Q96D46), NME2 (P22392), NME3
 (O13232), NME7 (Q9Y5B8), NMT1 (P30419), NNMT (P40261), NOB1 (Q9ULX3), NOC2L (Q9Y3T9), NOC3L
 (Q8WTT2), NOC4L (Q9BVI4), NOG (Q13253), NOL11 (Q9H8H0), NOL6 (Q9H6R4), NOL9 (Q5SY16), NOMO2
 (Q5JPE7), NONO (Q15233), NOP10 (Q9NPE3), NOP16 (Q9Y3C1), NOP2 (P46087), NOP56 (O00567), NOP58
 (Q9Y2X3), NOS1AP (O75052), NOSIP (Q9Y314), NOTCH2 (Q04721), NOVA2 (Q9UNW9), NPC1 (O15118),
 NPC2 (P61916), NPEPPS (P55786), NPLOC4 (Q8TAT6), NPM1 (P06748), NPTN (Q9Y639), NPW (Q8N729),
 NQO1 (P15559), NQO2 (P16083), NRAS (P01111), NRBP1 (Q9UHY1), NRD1 (O43847), NRP1 (O14786), NRP2
 (O60462), NSDHL (O15738), NSF (P46459), NSUN2 (Q08J23), NSUN5 (Q96P11), NSUN6 (Q8TEA1), NT5C
 (Q8TCD5), NT5C2 (P49902), NT5C3L (Q969T7), NT5E (P21589), NTN1 (O95631), NUBP1 (P53384), NUBP2
 (Q9Y5Y2), NUCB1 (Q02818), NUCKS1 (Q9H1E3), NUDC (Q9Y266), NUDCD1 (Q96RS6), NUDCD2 (Q8WVJ2),
 NUDT1 (P36639), NUDT10 (Q8NFP7), NUDT16 (Q96DE0), NUDT16L1 (Q9BRJ7), NUDT21 (O43809), NUDT4
 (Q9NZJ9), NUDT5 (Q9UKK9), NUMA1 (Q14980), NUP188 (Q5SRE5), NUP210 (Q8TEM1), NUP37 (Q8NFB4),
 NUP43 (Q8NFB3), NUP54 (Q7Z3B4), NUP62 (P37198), NUP85 (Q9BW27), NUP88 (Q99567), NUP93 (Q8N1F7),
 NUTF2 (P61970), NXF1 (Q9UBU9), NXN (Q6DKJ4), NXT1 (Q9UKK6), OAT (P04181), OBSL1 (O75147), OCRL
 (Q01968), ODR4 (Q5SWX8), ODZ2 (Q9NT68), ODZ3 (Q9P273), OGFOD1 (Q8N543), OGT (O15294), OLA1
 (Q9NTK5), OLFML3 (Q9NRN5), OPA1 (O60313), ORC3 (Q9UBD5), OSBP (P22059), OSBPL6 (Q9BZF3),
 OSGEP (Q9NPF4), OTUB1 (Q96FW1), OVCA2 (Q8WZ82), OXCT1 (P55809), OXSR1 (O95747), P4HA1
 (P13674), P4HB (P07237), PA2G4 (Q9UQ80), PAAF1 (Q9BRP4), PABPC1 (P11940), PABPC4 (Q13310),
 PABPN1 (Q86U42), PACSIN2 (Q9UNF0), PACSIN3 (Q9UKS6), PAF1 (Q8N7H5), PAFAH1B1 (P43034),
 PAFAH1B2 (P68402), PAFAH1B3 (Q15102), PAICS (P22234), PAIP1 (Q9H074), PAK1IP1 (Q9NWT1), PAK2
 (O13177), PALD (Q9ULE6), PALLD (Q8WX93), PANK4 (Q9NVE7), PAPOLA (P51003),

(continued)

PAPSS1 (Q43252), PARK7 (Q99497), PARN (Q95453), PARP1 (P09874), PARP4 (Q9UKK3), PARVA (Q9NVD7),
 PBLD (P30039), PCBD1 (P61457), PCBP1 (Q15365), PCBP2 (Q15366), PCDHB2 (Q9Y5E7), PCDHGC3
 (Q9UN70), PCID2 (Q5JVF3), PCMT1 (P22061), PCNA (P12004), PCOLCE2 (Q9UKZ9), PCYOX1 (Q9UHG3),
 PCYOX1L (Q8NBM8), PCYT2 (Q99447), PDCD10 (Q9BUL8), PDCD11 (Q14690), PDCD4 (Q53EL6), PDCD5
 (O14737), PDCD6 (O75340), PDCD6IP (Q8WUM4), PDCL3 (Q9H2J4), PDDC1 (Q8NB37), PDE12 (Q6L8Q7),
 PDGFRA (P16234), PDIA3 (P30101), PDIA4 (P13667), PDIA5 (Q14554), PDIA6 (Q15084), PDLIM1 (O00151),
 PDLIM4 (P50479), PDLIM5 (Q96HC4), PDLIM7 (Q9NR12), PDRO (Q6IAA8), PDS5A (Q29RF7), PDS5B
 (Q9NTI5), PDXK (O00764), PDXP (Q96GD0), PEA15 (Q15121), PEBP1 (P30086), PECI (O75521), PEF1
 (Q9UBV8), PELO (Q9BRX2), PELP1 (Q8IZL8), PEPD (P12955), PES1 (O00541), PFAS (O15067), PFDN1
 (O60925), PFDN2 (Q9UHV9), PFDN4 (Q9NQP4), PFDN5 (Q99471), PFDN6 (O15212), PFKL (P17858), PFKM
 (P08237), PFKP (Q01813), PFN1 (P07737), PFN2 (P35080), PGAM1 (P18669), PGAM5 (Q96HS1), PGD
 (P52209), PGGT1B (P53609), PGK1 (P00558), PGLS (O95336), PGLYRP2 (Q96PD5), PGM1 (P36871), PGM2L1
 (Q6PCE3), PGM3 (O95394), PGP (A6NDG6), PGRMC1 (O00264), PGRMC2 (O15173), PHB (P35232), PHB2
 (Q99623), PHF5A (Q7RTV0), PHF6 (Q8IWS0), PHGDH (O43175), PHKB (Q93100), PHLDA1 (Q8WV24),
 PHLDA3 (Q9Y5J5), PHLDB1 (Q86UU1), PHPT1 (Q9NRX4), PI15 (O43692), PI4KA (P42356), PICALM (Q13492),
 PIGT (Q969N2), PIK3CA (P42336), PIK3R4 (Q99570), PIN1 (O13526), PIP4K2A (P48426), PIP4K2B (P78356),
 PIP4K2C (Q8TBX8), PIPOX (Q9P0Z9), PIPSL (A2A3N6), PITPNB (P48739), PKM2 (P14618), PKP1 (O13835),
 PLAA (Q9Y263), PLCB3 (Q01970), PLCD1 (P51178), PLCD3 (Q8N3E9), PLCG1 (P19174), PLCG2 (P16885),
 PLD3 (Q8IV08), PLEC (Q15149), PLIN2 (Q99541), PLIN3 (O60664), PLK1 (P53350), PLOD1 (Q02809), PLOD2
 (O00469), PLOD3 (O60568), PLRG1 (O43660), PLS1 (O14651), PLS3 (P13797), PLSCR3 (Q9NRY6), PLTP
 (P55058), PLXNA1 (Q9UIW2), PLXNB2 (O15031), PLXND1 (Q9Y4D7), PMM2 (O15305), PMPCA (Q10713),
 PMPCB (O75439), PMVK (Q15126), PNMA2 (Q9UL42), PNN (Q9H307), PNO1 (Q9NRX1), PNP (P00491),
 PNPLA2 (Q96AD5), PODXL (O00592), POLD1 (P28340), POLD2 (P49005), POLE3 (Q9NRF9), POLR1A
 (Q95602), POLR1B (Q9H9Y6), POLR1C (O15160), POLR1D (Q9Y2S0), POLR2A (P24928), POLR2B (P30876),
 POLR2C (P19387), POLR2E (P19388), POLR2G (P62487), POLR2H (P52434), POLR2J (P52435), POLR2K
 (P53803), POLR3A (O14802), POLR3B (Q9NW08), POLR3C (Q9BUI4), POP1 (Q99575), POP4 (Q95707), POP7
 (O75817), POR (P16435), PPA1 (O15181), PPA2 (Q9H2U2), PPAN (Q9NQ55), PPAP2A (O14494), PPAT
 (Q06203), PPCS (Q9HAB8), PPFIBP1 (Q86W92), PPIA (P62937), PPIB (P23284), PPIC (P45877), PPID
 (Q08752), PPIF (P30405), PPIH (O43447), PPIL1 (Q9Y3C6), PPM1F (P49593), PPM1G (O15355), PPME1
 (Q9Y570), PPP1CA (P62136), PPP1CB (P62140), PPP1CC (P36873), PPP1R14B (Q96C90), PPP1R7 (Q15435),
 PPP1R8 (Q12972), PPP2CA (P67775), PPP2CB (P62714), PPP2R1A (P30153),

EP 3 033 418 B1

(continued)

PPP2R2A (P63151), PPP2R2D (Q66LE6), PPP2R4 (Q15257), PPP2R5D (Q14738), PPP2R5E (Q16537),
 PPP3CA (Q08209), PPP4C (P60510), PPP4R1 (Q8TF05), PPP5C (P53041), PPP6C (O00743), PPP6R3
 (Q5H9R7), PPPDE2 (Q6ICB0), PPT1 (P50897), PPWD1 (Q96BP3), PRCP (P42785), PRDX1 (Q06830), PRDX2
 (P32119), PRDX3 (P30048), PRDX4 (Q13162), PRDX6 (P30041), PREP (P48147), PREPL (Q4J6C6), PRIM1
 (P49642), PRIM2 (P49643), PRKAA1 (Q13131), PRKACA (P17612), PRKACB (P22694), PRKAG1 (P54619),
 PRKAR1A (P10644), PRKAR2A (P13861), PRKCA (P17252), PRKCI (P41743), PRKCSH (P14314), PRKDC
 (P78527), PRKRA (O75569), PRMT1 (Q99873), PRMT10 (Q6P2P2), PRMT3 (O60678), PRMT5 (O14744),
 PRMT7 (Q9NVM4), PROSC (O94903), PRPF19 (Q9UMS4), PRPF3 (O43395), PRPF31 (Q8WWY3), PRPF4
 (O43172), PRPF40A (O75400), PRPF4B (O13523), PRPF6 (O94906), PRPF8 (Q6P2Q9), PRPS1 (P60891),
 PRPS2 (P11908), PRPSAP2 (O60256), PRRC1 (Q96M27), PRSS23 (O95084), PRTFDC1 (Q9NRG1), PSAP
 (P07602), PSAT1 (Q9Y617), PSD3 (Q9NYI0), PSENEN (Q9NZ42), PSIP1 (O75475), PSMA1 (P25786), PSMA2
 (P25787), PSMA3 (P25788), PSMA4 (P25789), PSMA5 (P28066), PSMA6 (P60900), PSMA7 (O14818), PSMB1
 (P20618), PSMB2 (P49721), PSMB3 (P49720), PSMB4 (P28070), PSMB5 (P28074), PSMB6 (P28072), PSMB7
 (Q99436), PSMC1 (P62191), PSMC2 (P35998), PSMC3 (P17980), PSMC4 (P43686), PSMC5 (P62195), PSMC6
 (P62333), PSMD1 (Q99460), PSMD10 (O75832), PSMD11 (O00231), PSMD12 (O00232), PSMD13 (Q9UNM6),
 PSMD14 (O00487), PSMD2 (Q13200), PSMD3 (O43242), PSMD4 (P55036), PSMD5 (Q16401), PSMD6
 (Q15008), PSMD7 (P51665), PSMD8 (P48556), PSMD9 (O00233), PSME1 (Q06323), PSME2 (Q9UL46), PSME3
 (P61289), PSME4 (Q14997), PSMG1 (O95456), PSMG2 (Q969U7), PSPC1 (Q8WXF1), PSPH (P78330), PTBP1
 (P26599), PTGES2 (Q9H7Z7), PTGES3 (Q15185), PTGFRN (Q9P2B2), PTGR1 (Q14914), PTHLH (P12272),
 PTK2 (Q05397), PTK7 (Q13308), PTMA (P06454), PTN (P21246), PTP4A1 (Q93096), PTPN1 (P18031), PTPN11
 (Q06124), PTPN23 (Q9H3S7), PTPRA (P18433), PTPRE (P23469), PTPRG (P23470), PTPRJ (Q12913),
 PTPRZ1 (P23471), PUF60 (Q9UHX1), PURA (Q00577), PURB (Q96QR8), PUS1 (Q9Y606), PUS7 (Q96PZ0),
 PVR (P15151), PVRL2 (Q92692), PWP1 (Q13610), PWP2 (Q15269), PXDN (Q92626), P XK (Q7Z7A4), PXN
 (P49023), PYCR1 (P32322), PYCRL (Q53H96), PYGB (P11216), PYGL (P06737), OARS (P47897), QDPR
 (P09417), QKI (Q96PU8), QTRT1 (Q9BXR0), RAB10 (P61026), RAB11A (P62491), RAB11FIP1 (Q6WKZ4),
 RAB12 (Q6IQ22), RAB13 (P51153), RAB14 (P61106), RAB18 (Q9NP72), RAB1A (P62820), RAB1B (Q9H0U4),
 RAB21 (Q9UL25), RAB22A (Q9UL26), RAB23 (Q9ULC3), RAB27A (P51159), RAB2A (P61019), RAB2B
 (Q8WUD1), RAB32 (Q13637), RAB34 (Q9BZG1), RAB35 (Q15286), RAB3A (P20336), RAB3GAP1 (Q15042),
 RAB3GAP2 (Q9H2M9), RAB4A (P20338), RAB5A (P20339), RAB5B (P61020), RAB5C (P51148), RAB6A
 (P20340), RAB7A (P51149), RAB8A (P61006), RAB8B (Q92930), RABAC1 (Q9UI14), RABGAP1 (Q9Y3P9),
 RABGGTA (Q92696), RABGGTB (P53611), RABL2A (Q9UBK7), RABL3 (Q5HYI8), RAC1 (P63000), RAC3
 (P60763), RAD23B (P54727),

EP 3 033 418 B1

(continued)

RAD50 (Q92878), RAE1 (P78406), RAF1 (P04049), RALA (P11233), RALB (P11234), RALY (Q9UKM9), RAN (P62826), RANBP1 (P43487), RANBP2 (P49792), RANGAP1 (P46060), RAP1A (P62834), RAP1B (P61224), RAP1GDS1 (P52306), RAP2B (P61225), RAPH1 (Q70E73), RARS (P54136), RASA1 (P20936), RASA3 (Q14644), RBBP4 (Q09028), RBBP5 (Q15291), RBBP7 (Q16576), RBM12 (Q9NTZ6), RBM14 (Q96PK6), RBM15 (Q96T37), RBM22 (Q9NW64), RBM25 (P49756), RBM26 (Q5T8P6), RBM28 (Q9NW13), RBM39 (Q14498), RBM4 (Q9BWF3), RBM8A (Q9Y5S9), RBMX (P38159), RBP1 (P09455), RBPJ (Q06330), RBX1 (P62877), RCC1 (P18754), RCC2 (Q9P258), RCL (O43598), RCL1 (Q9Y2P8), RCN1 (Q15293), RDH11 (Q8TC12), RDH13 (Q8NBN7), RDX (P35241), RECQL (P46063), RELA (Q04206), REPS1 (Q96D71), RETSAT (Q6NUM9), RFC2 (P35250), RFC3 (P40938), RFC4 (P35249), RFC5 (P40937), RFFL (Q8WZ73), RFTN1 (Q14699), RHEB (Q15382), RHOA (P61586), RHOB (P62745), RHOC (P08134), RHOF (Q9HBH0), RHOG (P84095), RHOT2 (Q8IX11), RIC8A (Q9NPQ8), RNASEH2C (Q8TDP1), RNF114 (Q9Y508), RNF20 (Q5VTR2), RNF213 (Q63HN8), RNF7 (Q9UBF6), RNGTT (O60942), RNH1 (P13489), RNMT (O43148), RNPEP (Q9H4A4), ROBLD3 (Q9Y2Q5), ROCK1 (Q13464), ROCK2 (O75116), RP2 (O75695), RPA1 (P27694), RPA2 (P15927), RPA3 (P35244), RPE (Q96AT9), RPF2 (Q9H7B2), RPIA (P49247), RPL10 (P27635), RPL10A (P62906), RPL11 (P62913), RPL12 (P30050), RPL13 (P26373), RPL13A (P40429), RPL14 (P50914), RPL15 (P61313), RPL17 (P18621), RPL18 (Q07020), RPL18A (Q02543), RPL19 (P84098), RPL21 (P46778), RPL22 (P35268), RPL22L1 (Q6P5R6), RPL23 (P62829), RPL23A (P62750), RPL24 (P83731), RPL26 (P61254), RPL26L1 (Q9UNX3), RPL27 (P61353), RPL27A (P46776), RPL28 (P46779), RPL29 (P47914), RPL3 (P39023), RPL30 (P62888), RPL31 (P62899), RPL32 (P62910), RPL34 (P49207), RPL35 (P42766), RPL35A (P18077), RPL36 (Q9Y3U8), RPL36A (P83881), RPL36AL (Q969Q0), RPL37 (P61927), RPL37A (P61513), RPL38 (P63173), RPL4 (P36578), RPL5 (P46777), RPL6 (Q02878), RPL7 (P18124), RPL7A (P62424), RPL7L1 (Q6DKI1), RPL8 (P62917), RPL9 (P32969), RPLP0 (P05388), RPLP1 (P05386), RPLP2 (P05387), RPN1 (P04843), RPN2 (P04844), RPP30 (P78346), RPP38 (P78345), RPRD1A (Q96P16), RPRD1B (Q9NQG5), RPS10 (P46783), RPS11 (P62280), RPS12 (P25398), RPS13 (P62277), RPS14 (P62263), RPS15 (P62841), RPS15A (P62244), RPS16 (P62249), RPS17 (P08708), RPS18 (P62269), RPS19 (P39019), RPS2 (P15880), RPS20 (P60866), RPS21 (P63220), RPS23 (P62266), RPS24 (P62847), RPS25 (P62851), RPS26 (P62854), RPS27 (P42677), RPS27A (P62979), RPS27L (Q71UM5), RPS28 (P62857), RPS29 (P62273), RPS3 (P23396), RPS3A (P61247), RPS4X (P62701), RPS4Y1 (P22090), RPS5 (P46782), RPS6 (P62753), RPS6KA1 (Q15418), RPS6KA3 (P51812), RPS7 (P62081), RPS8 (P62241), RPS9 (P46781), RPSA (P08865), RQCD1 (Q92600), RRAGC (Q9HB90), RRAS2 (P62070), RRBP1 (Q9P2E9), RRM1 (P23921), RRM2 (P31350), RRM2B (Q7LG56), RRP1 (P56182), RRP12 (Q5JTH9), RRP1B (Q14684), RRP7A (Q9Y3A4), RRP9 (O43818), RRS1 (Q15050), RSL1D1 (O76021), RSL24D1 (Q9UHA3),

EP 3 033 418 B1

(continued)

RSPRY1 (Q96DX4), RSU1 (Q15404), RTCD1 (O00442), RTKN (Q9BST9), RTN3 (Q95197), RTN4 (Q9NQC3), RUVBL1 (Q9Y265), RUVBL2 (Q9Y230), RWDD2B (P57060), S100A10 (P60903), S100A11 (P31949), S100A13 (Q99584), S100A16 (Q96FQ6), S100A2 (P29034), S100A4 (P26447), S100A6 (P06703), S100A7 (P31151), S100A8 (P05109), S100A9 (P06702), SAAL1 (Q96ER3), SACS (Q9NZJ4), SAE1 (Q9UBE0), SAMHD1 (Q9Y3Z3), SAP18 (O00422), SAR1A (Q9NR31), SARM1 (Q6SZW1), SARNP (P82979), SARS (P49591), SARS2 (Q9NP81), SART3 (Q15020), SBDS (Q9Y3A5), SBF1 (O95248), SCARB1 (Q8WTV0), SCARB2 (Q14108), SCCPDH (Q8NBX0), SCFD1 (Q8WVM8), SCFD2 (Q8WU76), SCP2 (P22307), SCPEP1 (Q9HB40), SCRG1 (O75711), SCRIB (Q14160), SCRNI1 (Q12765), SCRNI2 (Q96FV2), SCYL1 (Q96KG9), SDC2 (P34741), SDC4 (P31431), SDCBP (O00560), SDCCAG1 (O60524), SDCCAG3 (Q96C92), SDHA (P31040), SDHB (P21912), SDK1 (Q7Z5N4), SDSL (Q96GA7), SEC13 (P55735), SEC14L2 (O76054), SEC22B (O75396), SEC23A (O15436), SEC23B (Q15437), SEC23IP (Q9Y6Y8), SEC24A (O95486), SEC24B (O95487), SEC24C (P53992), SEC24D (O94855), SEC31A (O94979), SEC61B (P60468), SEC61G (P60059), SEH1L (Q96EE3), SELH (Q8IZQ5), SELO (Q9BVL4), SEMA3A (Q14563), SENP3 (Q9H4L4), SEPSECS (Q9HD40), 40422 (Q9P0V9), 40787 (Q9NVA2), 37500 (Q15019), 38596 (Q99719), 39326 (Q16181), 40057 (Q9UHD8), SERBP1 (Q8NC51), SERPINB12 (Q96P63), SERPINB3 (P29508), SERPINB6 (P35237), SERPINH1 (P50454), SESN2 (P58004), SET (Q01105), SETD3 (Q86TU7), SF3A1 (Q15459), SF3A2 (Q15428), SF3A3 (Q12874), SF3B1 (O75533), SF3B14 (Q9Y3B4), SF3B2 (Q13435), SF3B3 (Q15393), SF3B4 (Q15427), SF3B5 (Q9BWJ5), SFN (P31947), SFPQ (P23246), SFRP4 (Q6FHJ7), SFXN3 (Q9BWM7), SGTA (O43765), SH3BGRL3 (Q9H299), SH3BP4 (Q9P0V3), SH3GL1 (Q99961), SH3GLB1 (Q9Y371), SHC1 (P29353), SHMT1 (P34896), SHMT2 (P34897), SHOC2 (Q9UQ13), SHPK (Q9UHH6), SIRT5 (Q9NXA8), SKIV2L (Q15477), SKIV2L2 (P42285), SKP1 (P63208), SLC12A2 (P55011), SLC12A4 (Q9UP95), SLC16A1 (P53985), SLC1A3 (P43003), SLC1A5 (Q15758), SLC25A10 (Q9UBX3), SLC25A11 (Q02978), SLC25A13 (Q9UJS0), SLC25A22 (Q9H936), SLC25A3 (Q00325), SLC25A5 (P05141), SLC25A6 (P12236), SLC26A2 (P50443), SLC29A1 (Q99808), SLC29A2 (Q14542), SLC2A1 (P11166), SLC30A1 (Q9Y6M5), SLC38A1 (Q9H2H9), SLC3A2 (P08195), SLC44A2 (Q8IWA5), SLC4A2 (P04920), SLC4A7 (Q9Y6M7), SLC5A3 (P53794), SLC5A6 (Q9Y289), SLC6A8 (P48029), SLC7A1 (P30825), SLC7A5 (Q01650), SLC9A3R1 (O14745), SLC9A3R2 (Q15599), SLIRP (Q9GZT3), SLK (Q9H2G2), SMAD1 (Q15797), SMAD2 (Q15796), SMARCA4 (P51532), SMARCA5 (O60264), SMARCB1 (Q12824), SMARCC1 (Q92922), SMARCC2 (Q8TAQ2), SMARCD2 (Q92925), SMC1A (Q14683), SMC2 (O95347), SMC3 (Q9UQE7), SMC4 (Q9NTJ3), SMC5 (Q8IY18), SMCHD1 (A6NHR9), SMEK1 (Q6IN85), SMG1 (Q96Q15), SMN1 (Q16637), SMS (P52788), SMU1 (Q2TAY7), SMYD3 (Q9H7B4), SMYD5 (Q6GMV2), SNAP23 (O00161), SND1 (Q7KZF4), SNF8 (Q96H20), SNRNP200 (O75643), SNRNP40 (Q96DI7), SNRNP70 (P08621), SNRPA1 (P09661), SNRPB (P14678),

EP 3 033 418 B1

(continued)

SNRPB2 (P08579), SNRPD1 (P62314), SNRPD2 (P62316), SNRPD3 (P62318), SNRPE (P62304), SNRPF (P62306), SNRPG (P62308), SNTB1 (Q13884), SNTB2 (Q13425), SNX1 (Q13596), SNX12 (Q9UMY4), SNX17 (Q15036), SNX18 (Q96RF0), SNX2 (O60749), SNX27 (Q96L92), SNX3 (O60493), SNX5 (Q9Y5X3), SNX6 (Q9UNH7), SNX9 (Q9Y5X1), SOD1 (P00441), SOD2 (P04179), SORD (Q00796), SORT1 (Q99523), SPATS2L (Q9NUQ6), SPC24 (Q8NBT2), SPCS2 (Q15005), SPCS3 (P61009), SPG21 (Q9NZD8), SPIN1 (Q9Y657), SPR (P35270), SPRR1B (P22528), SPRR2E (P22531), SPTAN1 (Q13813), SPTBN1 (Q01082), SPTBN2 (O15020), SR140 (O15042), SRBD1 (Q8N5C6), SRCRL (A1L4H1), SRGAP2 (O75044), SRI (P30626), SRM (P19623), SRP14 (P37108), SRP19 (P09132), SRP54 (P61011), SRP68 (Q9UHB9), SRP72 (O76094), SRP9 (P49458), SRPK1 (Q96SB4), SRPR (P08240), SRPRB (Q9Y5M8), SRPX (P78539), SRPX2 (O60687), SRR (Q9GZT4), SRRM1 (Q8IYB3), SRRM2 (Q9UQ35), SRRT (Q9BXP5), SRSF1 (Q07955), SRSF10 (O75494), SRSF11 (Q05519), SRSF2 (Q01130), SRSF3 (P84103), SRSF5 (Q13243), SRSF6 (Q13247), SRSF7 (Q16629), SRSF9 (Q13242), SRXN1 (Q9BYN0), SSB (P05455), SSBP1 (Q04837), SSR1 (P43307), SSR3 (Q9UNL2), SSRP1 (Q08945), SSSCA1 (O60232), SSU72 (Q9NP77), ST13 (P50502), STAG1 (Q8WVM7), STAM (Q92783), STAMPB (O95630), STAT1 (P42224), STAT2 (P52630), STAT3 (P40763), STAU1 (O95793), STIP1 (P31948), STK10 (O94804), STK24 (Q9Y6E0), STK25 (O00506), STK38 (Q15208), STK38L (Q9Y2H1), STOM (P27105), STOML2 (Q9UJZ1), STON2 (Q8WXE9), STRAP (Q9Y3F4), STT3A (P46977), STUB1 (Q9UNE7), STX12 (Q86Y82), STX4 (Q12846), STX5 (Q13190), STXBP1 (P61764), STXBP3 (O00186), STYX (Q8WUJ0), SUB1 (P53999), SUCLA2 (Q9P2R7), SUCLG2 (Q96I99), SUGT1 (Q9Y2Z0), SULF2 (Q8IWU5), SUMO1 (P63165), SUPT16H (Q9Y5B9), SUPT4H1 (P63272), SUPT5H (O00267), SUPT6H (Q7KZ85), SUSD5 (O60279), SVEP1 (Q4LDE5), SVIL (O95425), SWAP70 (Q9UH65), SYMPK (Q92797), SYNCRIP (O60506), SYNGR2 (O43760), SYNJ2BP (P57105), SYNM (O15061), SYPL1 (Q16563), TAB1 (Q15750), TAF9 (Q9Y3D8), TAGLN (Q01995), TAGLN2 (P37802), TALDO1 (P37837), TAOK1 (Q7L7X3), TARDBP (Q13148), TARS (P26639), TATDN1 (Q6P1N9), TAX1BP3 (O14907), TBC1D13 (Q9NVG8), TBC1D15 (Q8TC07), TBC1D23 (Q9NUY8), TBC1D24 (Q9ULP9), TBC1D4 (O60343), TBC1D9B (Q66K14), TBCA (O75347), TBCB (Q99426), TBCC (Q15814), TBCE (Q15813), TBK1 (Q9UHD2), TBL1XR1 (Q9BZK7), TBL2 (Q9Y4P3), TBL3 (Q12788), TBPL1 (P62380), TCEA1 (P23193), TCEB1 (Q15369), TCEB2 (Q15370), TCERG1 (O14776), TCF25 (Q9BQ70), TCP1 (P17987), TELO2 (Q9Y4R8), TEX10 (Q9NXF1), TEX15 (Q9BXT5), TF (P02787), TFCP2 (Q12800), TFG (Q92734), TERC (P02786), TGFB1 (P01137), TGFB2 (P61812), TGFB1 (Q15582), TGFB1 (Q15582), TGFB1 (Q15582), TGFBRAP1 (Q8WUH2), TGM1 (P22735), TGM3 (Q08188), TH1L (Q8IXH7), THBS1 (P07996), THBS3 (P49746), THG1L (Q9NWX6), THOC2 (Q8NI27), THOC3 (Q96J01), THOC5 (Q13769), THOC6 (Q86W42), THOC7 (Q6I9Y2), THOP1 (P52888), THTPA (Q9BU02), THUMPD1 (Q9NXG2), THUMPD3 (Q9BV44), THY1 (P04216), THYN1 (Q9P016),

(continued)

TIA1 (P31483), TIAL1 (Q01085), TIGAR (Q9NQ88), TIMM13 (Q9Y5L4), TIMM44 (O43615), TIMM50 (Q3ZCQ8),
 TIMM8A (O60220), TIMM8B (Q9Y5J9), TIMM9 (Q9Y5J7), TIMP2 (P16035), TIPRL (O75663), TJP1 (Q07157),
 5 TKT (P29401), TLN1 (Q9Y490), TLN2 (Q9Y4G6), TM9SF3 (Q9HD45), TMED10 (P49755), TMED2 (Q15363),
 TMED5 (Q9Y3A6), TMED7 (Q9Y3B3), TMED9 (Q9BVK6), TMEFF2 (Q9UIK5), TMEM132A (Q24JP5), TMEM2
 (Q9UHN6), TMEM30A (Q9NV96), TMEM33 (P57088), TMOD3 (Q9NYL9), TMPO (P42166), TMX1 (Q9H3N1),
 10 TNC (P24821), TNKS1BP1 (Q9C0C2), TNPO1 (Q92973), TNPO2 (O14787), TNPO3 (Q9Y5L0), TOM1L2
 (Q6ZVM7), TOMM20 (Q15388), TOMM34 (Q15785), TOMM5 (Q8N4H5), TOMM70A (O94826), TOP1 (P11387),
 TOP2A (P11388), TOP2B (Q02880), TP53I3 (Q53FA7), TP53RK (Q96S44), TPBG (Q13641), TPD52 (P55327),
 TPI1 (P60174), TPM1 (P09493), TPM2 (P07951), TPM3 (P06753), TPM3L (A6NL28), TPM4 (P67936), TPP2
 (P29144), TPT1 (P13693), TRA2A (Q13595), TRA2B (P62995), TRAF2 (Q12933), TRAP1 (Q12931), TRAPPC1
 (Q9Y5R8), TRAPPC2L (Q9UL33), TRAPPC3 (O43617), TRAPPC4 (Q9Y296), TRAPPC5 (Q8IUR0), TRIM16
 (O95361), TRIM22 (Q8IYM9), TRIM25 (Q14258), TRIM26 (Q12899), TRIM28 (Q13263), TRIM47 (Q96LD4),
 15 TRIM5 (Q9C035), TRIO (O75962), TRIP13 (Q15645), TRIP6 (Q15654), TRMT1 (Q9NXH9), TRMT112 (Q9UI30),
 TRMT5 (Q32P41), TRMT6 (Q9UJA5), TRMT61A (Q96FX7), TRNT1 (Q96Q11), TROVE2 (P10155), TRRAP
 (Q9Y4A5), TSG101 (Q99816), TSKU (Q8WUA8), TSN (Q15631), TSPAN14 (Q8NG11), TSPAN6 (O43657), TSR1
 (Q2NL82), TSSC1 (Q53HC9), TSTA3 (Q13630), TTC1 (Q99614), TTC15 (Q8WVT3), TTC27 (Q6P3X3), TTC37
 (Q6PGP7), TTC38 (Q5R3I4), TTC7B (Q86TV6), TTC9C (Q8N5M4), TTL (Q8NG68), TTLL12 (Q14166), TTN
 (Q8WZ42), TTYH1 (Q9H313), TTYH3 (Q9C0H2), TUBA1B (P68363), TUBA4A (P68366), TUBB (P07437),
 20 TUBB2B (Q9BVA1), TUBB2C (P68371), TUBB3 (Q13509), TUBB6 (Q9BUF5), TUBG1 (P23258), TUBGCP2
 (Q9BSJ2), TUBGCP3 (Q96CW5), TUFM (P49411), TWF1 (Q12792), TWF2 (Q6IBS0), TXN (P10599), TXNDC17
 (Q9BRA2), TXNDC5 (Q8NBS9), TXNDC9 (O14530), TXNL1 (O43396), TXNRD1 (Q16881), TYK2 (P29597),
 TYMS (P04818), U2AF1 (Q01081), U2AF2 (P26368), UAP1 (Q16222), UBA1 (P22314), UBA2 (Q9UBT2), UBA3
 (Q8TBC4), UBA52 (P62987), UBA6 (A0AVT1), UBE2D1 (P51668), UBE2D3 (P61077), UBE2E1 (P51965),
 25 UBE2G2 (P60604), UBE2I (P63279), UBE2J2 (Q8N2K1), UBE2K (P61086), UBE2L3 (P68036), UBE2M (P61081),
 UBE2N (P61088), UBE2O (Q9C0C9), UBE2S (Q16763), UBE2V1 (Q13404), UBE2V2 (Q15819), UBE3A
 (Q05086), UBE3C (Q15386), UBE4A (Q14139), UBE4B (O95155), UBF1 (O14562), UBL3 (O95164), UBL4A
 (P11441), UBL5 (Q9BZL1), UBLCP1 (Q8WVY7), UBP1 (Q9NZI7), UBQLN2 (Q9UHD9), UBR1 (Q8IUV7), UBR4
 (Q5T4S7), UBD1 (Q9HAC8), UBXN1 (Q04323), UBXN6 (Q9BZV1), UCHL1 (P09936), UCHL3 (P15374), UCHL5
 (Q9Y5K5), UCK2 (Q9BZX2), UFC1 (Q9Y3C8), UFD1L (Q92890), UGDH (O60701), UGGT1 (Q9NYU2), UGP2
 (Q16851), ULK3 (Q6PHR2), UMPS (P11172), UNC119B (A6NIH7), UNC45A (Q9H3U1), UPF1 (Q92900), UPP1
 (Q16831), UQCRC1 (P31930), UQCRC2 (P22695), UQCRFS1 (P47985),
 30 URB1 (O60287), URB2 (Q14146), UROD (P06132), UROS (P10746), USO1 (O60763), USP10 (Q14694), USP11
 (P51784), USP13 (Q92995), USP14 (P54578), USP15 (Q9Y4E8), USP24 (Q9UPU5), USP39 (Q53GS9), USP5
 (P45974), USP7 (Q93009), USP9X (Q93008), UTP15 (Q8TED0), UTP18 (Q9Y5J1), UTP20 (O75691), UTP6
 (Q9NYH9), UTRN (P46939), UXS1 (Q8NBZ7), UXT (Q9UBK9), VAC14 (Q08AM6), VAMP3 (Q15836), VAMP5
 (O95183), VAPA (Q9P0L0), VAPB (O95292), VARS (P26640), VASP (P50552), VAT1 (Q99536), VAV2 (P52735),
 VBP1 (P61758), VCAN (P13611), VCL (P18206), VCP (P55072), VDAC1 (P21796), VDAC2 (P45880), VDAC3
 (Q9Y277), VIM (P08670), VPRBP (Q9Y4B6), VPS11 (Q9H270), VPS13A (Q96RL7), VPS13C (Q709C8), VPS16
 (Q9H269), VPS18 (Q9P253), VPS24 (Q9Y3E7), VPS25 (Q9BRG1), VPS26A (O75436), VPS26B (Q4G0F5),
 40 VPS28 (Q9UK41), VPS29 (Q9UBQ0), VPS33A (Q96AX1), VPS33B (Q9H267), VPS35 (Q96QK1), VPS36
 (Q86VN1), VPS37B (Q9H9H4), VPS39 (Q96JC1), VPS41 (P49754), VPS45 (Q9NRW7), VPS4A (Q9UN37),
 VPS4B (O75351), VPS53 (Q5VIR6), VPS8 (Q8N3P4), VRK1 (Q99986), VTA1 (Q9NP79), VWA1 (Q6PCB0),
 VWA5A (O00534), WARS (P23381), WASF2 (Q9Y6W5), WASL (O00401), WBSCR22 (O43709), WDFY1
 (Q8IWB7), WDR1 (O75083), WDR11 (Q9BZH6), WDR12 (Q9GZL7), WDR18 (Q9BV38), WDR26 (Q9H7D7),
 WDR3 (Q9UNX4), WDR36 (Q8NI36), WDR4 (P57081), WDR43 (Q15061), WDR45L (Q5MNZ6), WDR48
 (Q8TAF3), WDR5 (P61964), WDR54 (Q9H977), WDR6 (Q9NNW5), WDR61 (Q9GZS3), WDR73 (Q6P4I2),
 WDR74 (Q6RFH5), WDR75 (Q8IWA0), WDR77 (Q9BQA1), WDR82 (Q6UXN9), WDR92 (Q96MX6), WHSC2
 (Q9H3P2), WRNIP1 (Q96S55), XP32 (Q5T750), XPC (Q01831), XPNPEP1 (Q9NQW7), XPO1 (O14980), XPO4
 (Q9C0E2), XPO5 (Q9HAV4), XPO6 (Q96QU8), XPO7 (Q9UIA9), XPOT (O43592), XRCC1 (P18887), XRCC5
 (P13010), XRCC6 (P12956), XRN2 (Q9H0D6), YARS (P54577), YBX1 (P67809), YES1 (P07947), YKT6 (O15498),
 45 YRDC (Q86U90), YTHDC1 (Q96MU7), YTHDF2 (Q9Y5A9), YWHAB (P31946), YWHAE (P62258), YWHAG
 (P61981), YWHAH (Q04917), YWHAQ (P27348), YWHAZ (P63104), ZC3H15 (Q8WU90), ZC3HAV1 (Q7Z2W4),
 ZC3HAV1L (Q96H79), ZCHC3 (Q9NUD5), ZFAND1 (Q8TCF1), ZFR (Q96KR1), ZMAT2 (Q96NC0), ZNF259
 (O75312), ZNF326 (Q5BKZ1), ZNF330 (Q9Y3S2), ZNF622 (Q969S3), ZNF765 (Q7L2R6), ZNFX1 (Q9P2E3),
 50 ZW10 (O43264), ZWILCH (Q9H900), ZYG11B (Q9C0D3), ZYX (Q15942).

Table 22: 100 most abundant proteins (name and SwissProt accession number) in CTX0E03 microvesicles

Identified proteins	Accession number
Actin, cytoplasmic 2	P63261
Histone H4	P62805
Histone H2B	Q99879
Histone H3.2	Q71DI3
Histone H2B type 1	P23527
Glyceraldehyde-3-phosphate dehydrogenase	P04406
Histone H2A type 2-A	Q6FI13
Ubiquitin-40S ribosomal protein S27a	P62979
Annexin A2	P07355
Alpha-enolase	P06733
Pyruvate kinase isozymes M1/M2	P14618
60S ribosomal protein L6	Q02878
Histone H2B type 2-E	Q16778
Heat shock cognate 71 kDa protein	P11142
Actin, alpha cardiac muscle 1	P68032
Heat shock protein HSP 90-beta	P08238
Histone H2B type 1-J	P06899
Elongation factor 1-alpha 1	P68104
Tubulin beta-2C chain	P68371
60S ribosomal protein L18	Q07020
Tubulin beta chain	P07437
40S ribosomal protein S2	P15880
40S ribosomal protein S11	P62280
Histone H2B type 3-B	Q8N257
Tubulin alpha-1 B chain	P68363
40S ribosomal protein S3	P23396
40S ribosomal protein S3a	P61247
Histone H2A type 1-D	P20671
Elongation factor 2	P13639
Heat shock protein HSP 90-alpha	P07900
GTP-binding nuclear protein Ran	P62826
60S ribosomal protein L4	P36578
40S ribosomal protein S9	P46781
Profilin-1	P07737
60S ribosomal protein L13a	P40429
Phosphoglycerate kinase 1	P00558
Fatty acid synthase	P49327

EP 3 033 418 B1

(continued)

	Identified proteins	Accession number
5	Annexin A1	P04083
	Histone H2A.Z	P0C0S5
	Vimentin	P08670
	40S ribosomal protein S6	P62753
10	Moesin	P26038
	Peptidyl-prolyl cis-trans isomerase A	P62937
	60S ribosomal protein L26	P61254
15	60S ribosomal protein L3	P39023
	40S ribosomal protein S8	P62241
	60S ribosomal protein L28	P46779
	Ezrin	P15311
20	40S ribosomal protein S4, X isoform	P62701
	60S ribosomal protein L7a	P62424
	60S ribosomal protein L13	P26373
25	60S ribosomal protein L7	P18124
	40S ribosomal protein S23	P62266
	60S ribosomal protein L5	P46777
	Eukaryotic initiation factor 4A-I	P60842
30	40S ribosomal protein S24	P62847
	Tubulin beta-2B chain	Q9BVA1
	60S ribosomal protein L8	P62917
35	60S ribosomal protein L15	P61313
	60S ribosomal protein L10	P27635
	Peroxiredoxin-1	Q06830
	Keratin, type I cytoskeletal 14	P02533
40	14-3-3 protein theta	P27348
	40S ribosomal protein S18	P62269
	Transketolase	P29401
45	60S ribosomal protein L24	P83731
	Histone H1.5	P16401
	Cofilin-1	P23528
	Dihydropyrimidinase-related protein 3	Q14195
50	60S ribosomal protein L21	P46778
	60S ribosomal protein L36	Q9Y3U8
	Sodium/potassium-transporting ATPase subunit alpha-1	P05023
55	40S ribosomal protein S16	P62249
	T-complex protein 1 subunit gamma	P49368
	Heterogeneous nuclear ribonucleoprotein A1	P09651

(continued)

Identified proteins	Accession number
60S ribosomal protein L14	P50914
Heat shock 70 kDa protein 1A/1B	P08107
T-complex protein 1 subunit theta	P50990
60S ribosomal protein L30	P62888
Protein S100-A6	P06703
40S ribosomal protein SA	P08865
CD44 antigen	P16070
60S ribosomal protein L35a	P18077
Tubulin beta-3 chain	Q13509
T-complex protein 1 subunit delta	P50991
4F2 cell-surface antigen heavy chain	P08195
T-complex protein 1 subunit beta	P78371
Myosin-9	P35579
Adenosylhomocysteinase	P23526
Filamin-A	P21333
Fatty acid-binding protein, brain	015540
Myristoylated alanine-rich C-kinase substrate	P29966
T-complex protein 1 subunit eta	Q99832
Fascin	Q16658
Fructose-bisphosphate aldolase A	P04075
60S ribosomal protein L27	P61353
60S ribosomal protein L17	P18621
Heterogeneous nuclear ribonucleoproteins A2/B1	P22626
60S ribosomal protein L10a	P62906
60S ribosomal protein L35	P42766

Discussion of proteomic data

[0359] CD63 (also known as MLA1 and TSPAN30), TSG101 (also known as ESCRT-I complex subunit TSG101), CD109 (also known as 150 kDa TGF-beta-1-binding protein) and thy-1 (also known as CD90) were detected in both exosomes and microvesicles.

[0360] Other tetraspanins were also detected: Tetraspanin-4, -5, -6, -9 and 14 were detected in the exosome fraction; tetraspanins-6 and -14 were detected in the microvesicles.

[0361] CD133 (also known as AC133, Prominin-1, PROM1, PROML1 and MSTP061) was detected in the exosomes but not the microvesicles.

[0362] CD53 (also known as MOX44 and TSPAN25), CD82 (also known as KAI1, SAR2, ST6 and TSPAN27), CD37 (also known as TSPAN26) and CD40 ligand (also known as CD40LG, CD40L and TNFSF5) were not detected in the exosomes or the microvesicles.

[0363] Nestin, GFAP and tubulin beta-3 chain (also known as TUBB3) were detected in both the exosome and microvesicle fractions, with tubulin beta-3 chain being particularly prominent within the top 100 proteins in both fractions. Sox2, DCX, GALC, GDNF and IDO were not detected.

[0364] Selectins and TNFR1 (also known as TNF receptor 1, TNFRSF1A, TNFAR and TNFR1) were not detected.

[0365] Integrin alpha-2, -3, -4, -5, -6, -7, -V and integrin beta-1, -4 and -8 were detected in both exosome and microvesicle

fractions. Integrin beta-3 and -5 were detected in the microvesicles only.

[0366] MHC Class I antigens (e.g. HLA_A1, HLA-A2 and HLA-B27) were detected in both the exosomes and microvesicles.

[0367] Cell-adhesion molecules (e.g. CADM1, CADM4, ICAM1, JAM3, L1CAM, NCAM) were detected in both the exosomes and microvesicles.

[0368] Cytoskeletal proteins (e.g. actin, vimentin, keratins, catenins, dystroglycan, neurofilament polypeptide, microtubule-associated protein, tubulin, desmoplaktin, plectin, plakophilin, septin, spectrin, talin, vinculin and zyxin) were detected in both the exosome and microvesicle fractions.

[0369] GTPases, clathrin, chaperones, heat-shock proteins (e.g. Hsp90, Hsp70), splicing factors, translation factors, annexins and growth factors (e.g. TGF-beta) were detected in both the exosomes and microvesicles.

[0370] Galectin-3, TIMP-1, thrombospondin-1, EGF receptor and CSK were detected in both the exosomes and microvesicles.

[0371] Figure 17 compares the proteomic data from the exosomes and microvesicles. Figure 17A illustrates the number of unique proteins within each micro particle population, isolated from week 2 Integra culture system. Figure 17B compares the biological processes associated with the identified proteins within each micro particle population, isolated from week 2 Integra system. The proteins identified within exosomes and microvesicles are associated with very similar biological processes.

[0372] Proteins associated with biotin metabolism were only found in exosomes and proteins involved in tryptophan biosynthesis and taurine/alpha-linolenic acid metabolism were only identified in microvesicles.

[0373] Figure 17C compares the CTX0E03 proteome to the Mesenchymal Stem Cell exosome proteome disclosed in Lai et al 2012, in which a total of 857 proteins were identified in exosomes released from mesenchymal stem cells.

[0374] Figure 17D compares the biological processes associated with the identified proteins within the MSC derived exosomes (Lim 2012) with the neural stem cell derived exosomes of the invention. The three biological processes found to be associated with the MSC derived exosomes only are (in decreasing order of significance): Asthma; phenylalanine, tyrosine and tryptophan biosynthesis; and primary immunodeficiency. The thirty biological processes found to be associated only with the neural stem cell derived exosomes are shown in Figure 18; the most significant biological function identified relates to RNA polymerase.

[0375] A further comparison of the 197 biological processes shared by both MSC derived exosomes and NSC derived exosomes shows that NSC exosomes contain notably more processes involved in RNA degradation, the Ribosome and spliceosomes, when compared to MSC exosomes.

[0376] The above comparison indicates a number of significant differences between NSC derived exosomes and MSC derived exosomes (as characterised by Lim *et al* 2012). The 4 most significant biological differences identified as present in NSC exosomes compared to being very low/absent in those identified by the Lim's group, all involve proteins associated with the production, packaging, function and degradation of genetic material, i.e RNA polymerase, RNA degradation, Ribosome and spliceosomes.

Example 20: Functional analysis of individual miRNAs

Methods

MiRNA mimic transfection and evaluation of cell proliferation by Cyquant

[0377] Twenty four hours prior to transfection, glioma cells, U373 or U87, were seeded into a 96-well plate. MiRNA transfection optimization was performed using AllStars Negative Control siRNA AF 488 (Qiagen). MiRNA transfection efficiency was 100% when the following conditions were used. 20 nM of each miRNA mimics (Qiagen), hsa-miR-1246 (SEQ ID No. 21), hsa-miR-4492 (SEQ ID no. 34), hsa-miR-4532 (SEQ ID No. 23), and hsa-miR-4488 (SEQ ID No. 61) were combined with Lipofectamine® 2000 (Invitrogen) and transfection performed according to manufacturer's instructions.

Experiment 1 (U373MG; 2500 cells/well; 10%FBS)

[0378] 2500 U373MG cells were seeded per 96-well and cultured in DMEM glutamax/10%FBS for 24hrs, 48hrs and 72hrs post-transfection. Cell proliferation was measured by CyQUANT® Cell Proliferation Assay Kit (Invitrogen). Briefly, following removal of the culture medium, 200 µl of the CyQUANT® GR dye/cell-lysis buffer was added into each well of the 96-well plate and incubated for 15 min. Fluorescence intensity of each well was obtained using a GloMax™ 96 microplate (Promega) plate counter at excitation and emission wavelengths of 480 and 520 nm, respectively.

Experiment 2 (U373MG; 2500 cell/well; 2%FBS)

[0379] 2500 U373MG cells were seeded per 96-well and cultured in DMEM glutamax/2%FBS for 24hrs, 48hrs and 72hrs post-transfection. Cell proliferation was measured by CyQUANT® Cell Proliferation Assay Kit (Invitrogen). Briefly, following removal of the culture medium, 200 µl of the CyQUANT® GR dye/cell-lysis buffer was added into each well of the 96-well plate and incubated for 15 min. Fluorescence intensity of each well was obtained using a GloMax™ 96 microplate (Promega) plate counter at excitation and emission wavelengths of 480 and 520 nm, respectively.

Experiment 3 (U87; 9000 cells/well; basal)

[0380] 9000 U87 cells were seeded per 96 well and cultured in EMEM+2nM glutamine for 0hrs, 24hrs, 48hrs and 72hrs post-transfection. Cell proliferation was measured by CyQUANT® Cell Proliferation Assay Kit (Invitrogen). Briefly, following removal of the culture medium, 200 µl of the CyQUANT® GR dye/cell-lysis buffer was added into each well of the 96-well plate and incubated for 15 min. Fluorescence intensity of each well was obtained using a GloMax™ 96 Microplate (Promega) plate counter at excitation and emission wavelengths of 480 and 520 nm, respectively.

Results

[0381] Next generation sequence (NGS) analysis of miRNA contents in CTX0E03-derived exosomes revealed the presence of a set of top-ranked miRNAs, hsa-mir-1246, hsa-mir-4488, hsa-mir-4492, and hsa-mir-4532. To assess the functionality of these individual miRNAs in reducing glioma cell proliferation, each (mimic) miRNA was transfected into two cell line models of glioma: U373MG and U87.

[0382] The incidence of reduction of cell proliferation was dependent on the glioma model and cell culturing conditions, but each of the four miRNAs tested significantly reduced tumour cell proliferation in at least one of the models. hsa-mir-4492 and hsa-mir-4532 significantly reduced cell proliferation in each of the models tested. The results of Experiments 1 to 3 are shown in Figure 23A, B and C, respectively.

Example 21: Tolerability and pilot efficacy of exosomes in U-87MG human glioblastoma subcutaneous xenografts.

OBJECTIVE

[0383] To assess the tolerability and pilot efficacy of exosomes in U-87MG subcutaneous xenografts

METHODS

[0384]

ANIMALS

Number	30
No. of Groups & No per Group	5 groups, 5 mice/group
Species	<i>Mus musculus</i>
Strain	Athymic nude (Hsd:AthymicNude-Fxn1 ^{nu})
Age	5-7 weeks
Gender	Female
Body Weight	N/A
Animal ID	Transponder chip according to PRECOS Standard Operating Methods (SOMs)
Acclimatisation	≥ 1 week
Implantation site	Left flank
Anaesthesia	In accordance with PRECOS SOMs
Housing	According to PPL 70/7317 and PRECOS Standard Operating Procedures (SOPs)
Animal	Harlan UK

CELL LINES, *IN VITRO* EXPANSION

5	Cell Line Name	U-87MG
	Supplier and catalogue number	ECACC, 89081402
	Culture Medium	EMEM culture medium (Sigma, UK) containing 10% (v/v) heat inactivated foetal bovine serum (Hyclone, Thermo Scientific, UK)
10	Number of mice to be implanted	30
	Number of cells per mouse	8x10 ⁶ per mouse
	Matrigel/Cultrex/PBS/other	PBS supplemented with 0.1% glucose
15	Volume of diluent per mouse	0.1mL
	Number of batches	1

[0385] Cells will be harvested, washed in the culture medium described above and cells with viability of $\geq 90\%$ will be re-suspended for *in vivo* administration. Cells will be stored on ice for a minimum period of time (e.g. no longer than 30 minutes) prior to implantation.

IMPLANTATION

[0386] *TRANSPONDER IMPLANTATION*: Implanted at initiation (tumour implantation).

[0387] *TUMOUR IMPLANTATION*: 8x10⁶ viable cells in 100 μ l PBS + 0.1% glucose will be injected subcutaneously into the left flank of each mouse. A total of 30 mice will be implanted.

DATA CAPTURE

[0388] Body weight, dosing and any comments relating to clinical condition will be captured in real-time using the study management software, StudyDirector (StudyLog Systems Inc.). Data will be exported into Microsoft Excel and/or Graph-Pad Prism for subsequent data analysis and transformation.

[0389] Study specific data capture schedules will be created in Excel and completed by the study team. These data capture schedules will include study specific clinical observations; the recording of these observations and will be included in the final report and uploaded to the study folder at the end of the study.

[0390] *BODY WEIGHT*: Mice will be weighed x3 weekly during the dosing phase, weekly thereafter; clinical condition monitored daily for the duration of the study by an experienced technician.

[0391] *TUMOUR MONITORING (inc. BLI)*: Tumour will be measured 3 times a week and tumour volumes will be estimated using the formula $0.5(L \times W^2)$ by measuring the tumour in two dimensions using electronic callipers for the duration of the study.

[0392] *TREATMENT INITIATION AND DURATION*: The mice will be randomly allocated to the treatment groups (e.g. using a stratified randomisation software tool) such that there is a similar distribution of tumour size within and between treatment groups. Dosing will be initiated when the mean tumour volume of groups approximates 100-150mm³. The study will terminate 3 weeks following initiation.

TEST & REFERENCE SUBSTANCE ID STORAGE & FORMULATION

[0393]

Test & Reference Substance ID and storage.

Compound ID	Compound Source	Compound Storage	Vehicle Name	Vehicle Source	Vehicle Storage	Post-formulation storage
Exosome 0	ReNeuron Ltd	-80°C	0.9% saline	ReNeuron Ltd.	2-8°C	2-8°C Kept on ice during administration

EP 3 033 418 B1

(continued)

Compound ID	Compound Source	Compound Storage	Vehicle Name	Vehicle Source	Vehicle Storage	Post-formulation storage
Temozolomide	PRECOS Ltd	RT powder	10% DMSO	PRECOS	n/a	+4°C

[0394] Prepare dosing solutions freshly made before dosing.

DOSING

[0395] Mice will be dosed according to the following dosing schedule:

Group (No per group)	Compound ID	Dose mg/kg	Dose Volume (mg/ml)	Dose conc. (mg/ml)	Route	Dosing Frequency (bid/qd/tid) including wording e.g. twice daily etc.
1(5)	Vehicle	n/a	50µl ¹	0	Intratumoural	Once only
2(5)	Exosome 0	1	50µl ¹	TBD ²	Intratumoural	Once only
3(5)	Exosome 0	0.5	50µl ¹	TBD ²	Intratumoural	Once only
4(5)	Exosome 0	0.1	50µl ¹	TBD ²	Intratumoural	Once only
5(5)	Temozolomide	5	10.0	5	p.o.	Daily (q.d.)
¹ Dose calculated on mean body weight and delivered in a fixed volume of 50µl.						
² Concentration to be determined, dependent upon mean group body weight.						

STUDY ENDPOINTS/BODY WEIGHT LOSS (BWL) DURING THE STUDY

[0396]

- Terminate any mouse with sudden body weight loss approaching 20%
- Any mouse with continuous BWL approaching 20% over several daily measurements will be removed and terminated.

[0397] After one measurement of **body weight loss (BWL) > 10%**, a dose holiday will be given to the individual mouse. All dose holidays **must** be recorded on a Protocol Deviation.

[0398] *Whether to give dose holidays to all the mice or the individual mouse in the group should be done so in consultation with the client, but is ultimately at the Named Persons (or appointed deputies) discretion based on the severity/incidence of the BWL.*

TERMINATION

[0399] Each mouse will remain in the study until terminated (day 21), or until circumstances necessitate removal of an animal from the study e.g. loss of clinical condition and/or body weight.

[0400] Animals may also be terminated at any time during the study if any adverse effects are noted according to Home Office Project Licence PPL 70/7317.

[0401] Termination will be performed in accordance with United Kingdom Home Office Animals, (Scientific Procedures) Act 1986 and PRECOS SOPs.

TERMINAL SAMPLES

[0402] At termination the tumour will be excised and weighed. Tumours will be fixed in 10% Neutral Buffered Formalin and processing to FFPE blocks.

ANIMAL WELFARE AND REGULATION GUIDELINES*Housing And Environment*

[0403] Mice will be housed and cared for in accordance with the UK Animals (Scientific Procedures) Act 1986 (ASPA) and in line with the Directive 2010/63/EU of the European Parliament and of the Council of 22 September 2010 "on the protection of animals used for scientific purposes" and according to the and PRECOS Policies, SOPs and SOMs.

Animal Welfare Monitoring

[0404] This study will be conducted in line with the FELASSA Guidelines on Pain & Suffering in Experimental Animals and the NCRI Guidelines for the welfare and use of animals in Cancer Research (Workman et al., British Journal of Cancer (2010) 102, 1555-1577).

[0405] An experienced technician will check the condition of the mice at least daily. Unexpected adverse effects will be recorded and reported to the Named Animal Care & Welfare Officer (NACWO) and Named Veterinary Surgeon (NVS).

[0406] Animals may be terminated at any time during the study if any unexpected adverse effects are noted according to Home Office Project Licence PPL 70/7317 and the permitted severity band.

STATISTICS AND REPORTING*Statistical methods*

[0407] Statistical analysis if required will be performed in appropriate using the Minitab or PRISM statistical programmes for the PC.

RESULTS

[0408] Tumours were implanted on day 0 and measured from day 6; tumours were measured three times weekly by callipers and the tumour volume calculated. When tumours reached a mean tumour volume of ~165mm³ they were assigned to treatment groups based on mean tumour volume per cage in order to achieve a minimum amount of variation between and within groups.

[0409] The individual tumour volumes of each group on the day of assignment are presented in Figure 24; the mice were dosed on study day 12. The raw data for individual tumour volumes can be found in Figure 32.

[0410] Mouse body weights were monitored for the duration of the study. The data, expressed as the mean + standard error of the mean (% of the pre-dose weight), is presented graphically in Figure 25 (the dotted vertical line indicates the commencement of the dosing phase); the raw data for individual body weights can be found in Figures 30 and 31, absolute and relative measurements respectively. Body weight was stable over the duration of the study for each of the test agents and no adverse effects relating to the dosing protocols were documented in any of the treatment groups.

[0411] Mouse IDs recruited in dosing groups 1-4 received a single dose of increasing dose levels of Exosome 0 on study day 12. The Temozolomide dosing phase continued until study day 46 (35 oral doses). However, a number of mice were terminated prior to this point due to a number of listed adverse effects (see Figure 34; tumours reached the maximum permitted size as defined by UKCCCR guidelines (mean diameter 15mm)).

[0412] Figure 26 summarises the mean tumour volume for the treatment groups measured during the study, expressed as the group mean + standard error of the mean (% of the pre-dose volume).

[0413] One tumour from vehicle group 1 (ID4; 00077E7FDB-14) failed to demonstrate progressive growth following assignment and regressed to zero volume by day 29, as this is an untreated group the mouse has been classified as an outlier and removed from all analysis.

[0414] Loss of mice due to early terminations results in a reduction of the mean tumour volume from day 27 onwards. Figure 27 displays the tumour volume data (group mean + standard error of the mean; % pre-dose volume) of Figure 26 but in a truncated format *i.e.* all the line plots are graphed up to study day 25 before termination as a result of adverse effects occurred. The raw data for individual tumour volumes and individual tumour plots are detailed in appendices 3 and 4 respectively.

[0415] Mean tumour volumes were analysed statistically using a two-way ANOVA test (Figure 27 data set; GraphPad Prism; GraphPad Software, Inc.) for day 25. (The tumour volume data of G1 mouse ID14 was excluded from the statistical analysis; individual TV plots detailed in Appendix 4.) Although **there was a trend showing reduction in tumour volume for both 1mg/kg Exosome 0** (group 2) and Temozolomide (group 5), no statistically significant reduction in tumour volume was observed when compared to the vehicle group over the course of the study (GraphPad Prism; two-way ANOVA). The Bonferroni multiple comparison post-test did indicate a statistically significant reduction in tumour volume

on day 25 for group 5 versus group 1 ($p > 0.01$).

[0416] Terminal tumour weights were analysed statistically using a one-way ANOVA test (GraphPad Prism; GraphPad Software, Inc.); individual group comparisons were carried out on the total group tumour weights. No statistically significant differences among the mean tumour weight for each treatment group were observed using one-way ANOVA ($p = 0.2703$).

[0417] From the final tumour weight assessment (Figure 28; expressed as group mean + standard error of the mean (tumour weight)), **what is noticeable in the Exosome 0 1mg/kg dosing group (group 1) is that some of the tumours showed sensitivity to the treatment (IDs 12 and 21)**. Additionally, Temozolomide appears to increase the latency of the tumour instead of a significant decrease in tumour volume.

[0418] In survival analysis (Figure 29) utilising mean tumour diameter (15mm) as the humane survival endpoint, a **trend in increased survival was observed for 1mg/kg Exosome 0** and Temozolomide. However, no significant increase in survival for any of the treatment groups was observed when compared with the vehicle group ($p = 0.3651$; Log-Rank (Mantel-Cox) test; GraphPad Prism; GraphPad Software, Inc.). Temozolomide did result in an increase in survival compared with the two lower doses of Exosome 0 (groups 3 and 4; $p \geq 0.05$).

DISCUSSION

[0419] The primary objective of this pilot study was to assess the effect of several doses of Exosome 0 on the growth of U87MG subcutaneous glioblastoma xenografts and a dose of Temozolomide; an oral alkylating agent commonly used for the treatment of glioblastoma multiforme.

[0420] The agents under test in this study were well tolerated with no loss of body weight or adverse effects relating to treatment noted; however, tumour size and ulceration resulted in a decrease in mice per group from day 27. As the group sizes in this pilot study were already small the statistical significance that could be achieved with the test agents was therefore limited. The three dosing levels of the Exosome 0 were inefficacious in significantly reducing tumour volume (or tumour weight) when compared to the vehicle group, however, 40% of the tumours treated with the highest dose of Exosome 0 (1mg/kg; Group 1) showed sensitivity to the treatment.

[0421] Similarly the reduction in tumour size with Temozolomide was not significant, which also suggest the group sizes were too small to achieve statistical significance.

[0422] A trend in increased survival was also observed for 1mg/kg Exosome 0 and Temozolomide however significance was not achieved *versus* the vehicle group.

[0423] In conclusion, the dose levels of the Exosome 0 used in this study were well tolerated, but efficacy was emerging but not significant which was confounded by group size. Samples collected from this study could be analysed further for effect on proliferation, angiogenesis, necrosis and apoptosis. Further investigation using larger group sizes and higher dose of Exosome, if tolerated and soluble, could yield significant results.

Example 22: Histological Evaluation of Slides for U-87MG human glioblastoma subcutaneous xenografts

Summary

[0424] Tissues for histopathological examination (from Example 21) were stained with haematoxylin and eosin before being subjected to histopathological evaluation. This examination was to determine any differences in the appearance of U87 human glioblastoma tumours in animals given Exosome 0 or Temozolomide when compared to those given a vehicle alone.

[0425] In one animal given 1mg/kg Exosome 0 there was a particularly dramatic and effective ablation of the tumour mass.

Study Aims

[0426] The study was designed to investigate the properties of Exosome 0 in an in vivo model of tumorigenesis. This study was an investigation of the activity of this product in vivo, to assess tolerability and compare with an existing agent (temozolomide).

[0427] Temozolomide is an oral chemotherapy drug. It is an alkylating agent used for the treatment of glioblastoma. Temozolomide is also indicated for relapsed Grade III anaplastic astrocytoma, replacing the less well tolerated PCV (Procarbazine-Lomustine-Vincristine) regimen.

Methods

[0428] Histological Examination to determine any differences in the appearance of U87 human glioblastoma tumours in animals given Exosome 0 or Temozolomide when compared to those given a vehicle alone.

Results

1. Microscopic Findings

[0429] The tumours examined were large ovoid masses of apparently comparable sizes in the majority of cases, although the masses in animals 12 (Exosome 1.0 mg/kg), 1 and 9 (Temozolomide) were noticeably smaller. The majority of the tumours had necrotic centres and other necrotic foci within the mass. The extent of the necrosis was quite variable and did affect the appearance of the tumours, but there was no clear difference in the extent of necrosis between the groups. The tumour cells themselves were clonal with a little dysplasia and occasional apoptotic cells. The mitotic rate was relatively low and appeared to be consistent between treatment groups. The response of the host seemed variable, with several showing no evidence of a host response at all while others showed a slight to moderate inflammatory response with some fibroplasia in occasional animals forming a rudimentary capsule.

a) Vehicle Controls

[0430] The tumours in the vehicle control group all had the appearance as described above. The extent of central necrosis was minimal in animal 20, but in the rest of the animals was fairly extensive. The inflammatory response in 17 was greater than in the other members of this group.

b) Exosome 0 (1mg/kg)

[0431] The majority of tumours had an appearance that was indistinguishable from the tumours that were seen in the vehicle control animals. There was however one animal (12) where there was a dramatic response. In this animal the tumour appeared to have completely infarcted and there were no viable tumour cells visible in the section presented, only dense fibrous tissue and a slight infiltration of inflammatory cells, a large proportion of which appeared to be macrophages.

c) Exosome 0 (0.5 mg/kg)

[0432] The tumours in this group had an appearance that was indistinguishable from the tumours that were seen in the vehicle control animals.

d) Exosome 0 (0.1 mg/kg)

[0433] The tumours in this group had an appearance that was indistinguishable from the tumours that were seen in the vehicle control animals.

e) Temozolomide (5 mg/kg)

[0434] Three of the tumours in this group had an appearance that was indistinguishable from the tumours that were seen in the vehicle control animals, but there were two animals (1 and 9) where there seemed to be a response to treatment. In both animals the tumour had shrunk to quite a small size but there still appeared to be a substantial number of tumour cells in the section. These cells displayed rather more atypia than was seen in the other tumours, perhaps indicating an selective killing of the majority of cells, but potentially a selection of a more malignant phenotype by the test item. In two of the remaining animals (19 and 25) in this group tumours appeared similar in size to those seen in the vehicle control, there was however a clearly reduced cellularity compared to the vehicle control groups, but the increased atypia seen in animals 1 and 9 was again apparent. In the remaining animal the appearance was similar to the vehicle control.

2. Discussion and Conclusion

[0435] The consistent appearance of the tumour indicates a robust test system, which is suitable for assessment of efficacy. In one animal given 1mg/kg Exosome 0 there was a particularly dramatic and effective ablation of the tumour mass. In the animals that had received Temozolomide effects were seen in more animals, but the long term efficacy is perhaps more questionable as the effect appeared to be the selection of more atypical tumour cells, potentially with resistance to Temozolomide.

[0436] Given the consistency of the tumour appearance, it would suggest that the genotype of the tumours is well preserved and that, without being bound by theory - the large difference seen in animal 12 may be a result of a specific

Exosome/host interaction, rather than a direct effect of Exosome 0 on the tumour.

REFERENCES

[0437]

- Ambros et al RNA 2003. 9: 277-279
- Banerjee, S., Williamson, D., Habib, N., Gordon, M., Chataway, J. (2011) Age and Ageing 40:7-13
- Chung et al., Cell Stem Cell, 2, 113-117, 2008
- Dai, L.J., Moniri, M.R., Zeng, Z.R. et al. (2011) Cancer Lett 305(1):8-20.
- Ding, D. C., Shyu, W. C., Lin S. Z. (2011) Cell Transplant 20: 5-14
- Einstein, O., Ben-Hur, T. (2008) Arch Neurol 65:452-456
- Gennaro (2000) Remington: The Science and Practice of Pharmacy. 20th edition, ISBN: 0683306472
- Hassani Z, O'Reilly J, Pearse Y, Stroemer P, Tang E, Sinden J, Price J, Thuret S. "Human neural progenitor cell engraftment increases neurogenesis and microglial recruitment in the brain of rats with stroke." PLoS One. 2012;7(11):e50444. doi: 10.1371/journal.pone.0050444. Epub 2012 Nov 21.
- Hodges et al. Cell Transplant. 2007;16(2):101-15
- Horie, N., Pereira, N.P., Niizuma, K. Sun, G. et al. (2011) Stem Cells 29:274-285.
- Hulkower, K. I., Herber, R. L., (2011) Pharmaceuticals 3:107-124
- Katsuda, Kosaka, Takeshita, Ochiya. Proteomics 2013, 00, 1-17
- Katsuda, Tsuchiya, Kosaka, Yoshioka, Takagaki, Oki, Takeshita, Sakai, Kuroda, Ochiya. Scientific Reports 2013, 3:1197, p1-11.
- Kornblum, Stroke 2007, 38:810-816
- Lai et al "Proteolytic Potential of the MSC Exosome Proteome: Implications for an Exosome-Mediated Delivery of Therapeutic Proteasome". International Journal of Proteomics (2012) Article ID 971907, 14 pages.
- Littlewood, T. D., Hancock, D. C., Danielian, P. S. et al. (1995) Nucleic Acid Research 23:1686-1690.
- Miljan, E.A. & Sinden, J.D. (2009) Current Opinion in Molecular Therapeutics 4:394-403
- Miljan EA, Hines SJ, Pande P, Corteling RL, Hicks C, Zbarsky V, Umachandran, M, Sowinski P, Richardson S, Tang E, Wieruszew M, Patel S, Stroemer P, Sinden JD. Implantation of c-mycER TAM immortalized human mesencephalic-derived clonal cell lines ameliorates behavior dysfunction in a rat model of Parkinson's disease. Stem Cells Dev. 2009 Mar;18(2):307-19
- Mitchell et al Journal of Immunological Methods 335 (2008) 98-105
- Pollock et al, Exp Neurol. 2006 May;199(1):143-55.
- Mark F Pittenger; Alastair M Mackay; Stephen C Beck; Rama K Jaiswal; et al Science; Apr 2, 1999; 284, 5411
- Shah, K., (2012) Adv Drug Deliv Rev 64(8):739-748.
- Smith, E. J., Stroemer, R.P., Gorenkova, N., Nakajima, M. et al. (2012) Stem Cells 30:785-796.
- Stevenato, L., Corteling, R., Stroemer, P., Hope, A. et al. (2009) BMC Neuroscience 10:86
- Stroemer, P., Patel, S., Hope, A., Oliveira, C., Pollock, K., Sinden, J. (2009) Neurorehabil Neural Repair 23: 895-909.
- Théry, C., Ostrowski, M., Segura, E. et al. (2009) Nature Reviews Immunology 9: 581-593
- Their et al, "Direct Conversion of Fibroblasts into Stably Expandable Neural Stem Cells". Cell Stem Cell. 2012 Mar 20.
- Timmers, L., Lim, S. K., Arslan, F., Armstrong, J. S. et al. (2007) Stem Cell Res 1: 129-137
- Yuan, S.J., Martin, J., Elia, J., Flippin, J. et al. (2011) PLoS ONE 6:e17540

SEQUENCE LISTING

[0438]

<110> Reneuron Limited

<120> Stem Cell Microparticles and miRNA

<130> P062223WO

<141> 2014-02-12

<150> GB 1314573.5

<151> 2013-08-14

<150> GB 1317887.6

<151> 2013-10-09

<160> 752

<170> SeqWin2010, version 1.0

<210> 1

<211> 22

<212> RNA

<213> Homo sapiens

<400> 1

ugagguagua gguuguauag uu 22

<210> 2

<211> 23

<212> RNA

<213> Homo sapiens

<400> 2

uaccuguag auccgaauu gug 23

<210> 3

<211> 22

<212> RNA

<213> Homo sapiens

<400> 3

aaccguaga uccgaacuug ug 22

<210> 4

<211> 22

<212> RNA

<213> Homo sapiens

<400> 4

caccguaga accgaccuug cg 22

<210> 5

<211> 22

<212> RNA

<213> Homo sapiens

<400> 5

uccuguacug agcugccccg ag 22

<210> 6

<211> 21

<212> RNA

<213> Homo sapiens

<400> 6

uucacagugg cuaaguucug c 21

<210> 7

<211> 22

<212> RNA

	<213> Homo sapiens	
5	<400> 7 uauugcacuu gucccgccu gu	22
	<210> 8 <211> 23 <212> RNA <213> Homo sapiens	
10	<400> 8 caacggauc caaaagcag cug	23
15	<210> 9 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 9 uagcuauca gacugauguu ga	22
25	<210> 10 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 10 ugagguagua aguuguauug uu	22
	<210> 11 <211> 22 <212> RNA <213> Homo sapiens	
35	<400> 11 ugagguagua gauuguauag uu	22
40	<210> 12 <211> 22 <212> RNA <213> Homo sapiens	
45	<400> 12 uucaaguaau ccaggauagg cu	22
50	<210> 13 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 13 uauugcacuc gucccgccu cc	22
55	<210> 14 <211> 22 <212> RNA <213> Homo sapiens	

	<400> 14	
	ucggauccgu cugagcuugg cu	22
5	<210> 15	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 15	
	aacauucaac gcugucggug agu	23
	<210> 16	
	<211> 24	
15	<212> RNA	
	<213> Homo sapiens	
	<400> 16	
	uuuggcaaug guagaacuca cacu	24
20	<210> 17	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 17	
	ugagguagua gguuguaugg uu	22
	<210> 18	
30	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 18	
35	ugguagacua uggaacguag g	21
	<210> 19	
	<211> 22	
	<212> RNA	
40	<213> Homo sapiens	
	<400> 19	
	ugagaacuga auuccauagg cu	22
	<210> 20	
45	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 20	
50	caacaccagu cgaugggcug u	21
	<210> 21	
	<211> 19	
	<212> RNA	
55	<213> Homo sapiens	
	<400> 21	
	aauggauuuu uggagcagg	19

	<210> 22	
	<211> 22	
	<212> RNA	
5	<213> Homo sapiens	
	<400> 22	
	ugagguagua guuugugcug uu	22
10	<210> 23	
	<211> 17	
	<212> RNA	
	<213> Homo sapiens	
15	<400> 23	
	ccccggggag cccggcg	17
20	<210> 24	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 24	
	uauggcacug guagaauca cu	22
25	<210> 25	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
30	<400> 25	
	cuagacugaa gcuccuugag g	21
35	<210> 26	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
40	<400> 26	
	aaugcacccg ggcaaggauu cu	22
45	<210> 27	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 27	
	ugagguagga gguuguauag uu	22
50	<210> 28	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55	<400> 28	
	ugagguagua gguugugugg uu	22
	<210> 29	
	<211> 22	

	<212> RNA	
	<213> Homo sapiens	
5	<400> 29 uagcagcacg uaaauauugg cg	22
	<210> 30	
	<211> 22	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 30 uguaaacauc cucgacugga ag	22
15	<210> 31	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 31 uguaaacauc cccgacugga ag	22
	<210> 32	
	<211> 23	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 32 agguuacccg agcaacuug cau	23
30	<210> 33	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 33 aagcugccag uugaagaacu gu	22
	<210> 34	
40	<211> 17	
	<212> RNA	
	<213> Homo sapiens	
	<400> 34 ggggcugggc gcgcgcc	17
45	<210> 35	
	<211> 24	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 35 ucccugagac ccuuuaccu guga	24
55	<210> 36	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	

	<400> 36 agcuacaucu ggcucacuggg u	21
5	<210> 37 <211> 21 <212> RNA <213> Homo sapiens	
10	<400> 37 ucgaggagcu cacagucuag u	21
15	<210> 38 <211> 23 <212> RNA <213> Homo sapiens	
20	<400> 38 aacauucauu gcugucggug ggu	23
25	<210> 39 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 39 accuggcaua caauguagau uu	22
35	<210> 40 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 40 caaagaauuc uccuuuuggg cu	22
45	<210> 41 <211> 23 <212> RNA <213> Homo sapiens	
50	<400> 41 ugaggggcag agagcgagac uuu	23
55	<210> 42 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 42 ucccugagac ccuaacuugu ga	22
	<210> 43 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 43 ugagguagua guuuguacag uu	22

	<210> 44	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 44	
	augcaccugg gcaaggauuc ug	22
10		
	<210> 45	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 45	
	uguaaacauc cuugacugga ag	22
20		
	<210> 46	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 46	
	uucacagugg cuaaguuccg c	21
25		
	<210> 47	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 47	
	gaauguugcu cggugaaccc cu	22
35		
	<210> 48	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 48	
	ucagugcauc acagaacuuu gu	22
40		
	<210> 49	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
45		
	<400> 49	
	acggguuagg cucuugggag cu	22
50		
	<210> 50	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 50	
	aaauaacaac agauggccug u	21
	<210> 51	
	<211> 22	

	<212> RNA	
	<213> Homo sapiens	
5	<400> 51 uauacaaggg caagcucucu gu	22
	<210> 52	
	<211> 22	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 52 aaccgguaga uccgaucuug ug	22
15	<210> 53	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 53 agagguagua gguugcauag uu	22
	<210> 54	
	<211> 22	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 54 ucagugcacu acagaacuuu gu	22
30	<210> 55	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 55 aucacauugc cagggaauuc c	21
	<210> 56	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 56 cacuagauug ugagcuccug ga	22
	<210> 57	
	<211> 23	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 57 agcucggucu gaggccccuc agu	23
55	<210> 58	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	

	<400> 58	
	ucuuuggguua ucuagcugua uga	23
5	<210> 59	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 59	
	aucacauugc cagggauuac c	21
	<210> 60	
	<211> 23	
	<212> RNA	
15	<213> Homo sapiens	
	<400> 60	
	cacccggcug ugugcacaug ugc	23
20	<210> 61	
	<211> 18	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 61	
	agggggcggg cuccggcg	18
	<210> 62	
	<211> 23	
30	<212> RNA	
	<213> Homo sapiens	
	<400> 62	
	agcagcauug uacagggcua uga	23
35	<210> 63	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
40	<400> 63	
	cauugcacuu gucucggucu ga	22
	<210> 64	
45	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 64	
50	uuauaucgg acaaccuug u	21
	<210> 65	
	<211> 21	
	<212> RNA	
55	<213> Homo sapiens	
	<400> 65	
	acuggacuug gagucagaag g	21

<210> 66
 <211> 23
 <212> RNA
 <213> Homo sapiens
 5
 <400> 66
 uguaaacauc cuacacucuc age 23
 10
 <210> 67
 <211> 20
 <212> RNA
 <213> Homo sapiens
 15
 <400> 67
 uaacggccgc gguaccuuaa 20
 20
 <210> 68
 <211> 22
 <212> RNA
 <213> Homo sapiens
 25
 <400> 68
 ucacaaguca ggcucuuggg ac 22
 30
 <210> 69
 <211> 21
 <212> RNA
 <213> Homo sapiens
 35
 <400> 69
 uccgguucuc agggcuccac c 21
 40
 <210> 70
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 70
 uuaucagaau cuccaggggu ac 22
 50
 <210> 71
 <211> 22
 <212> RNA
 <213> Homo sapiens
 55
 <400> 71
 cuuucagucg gauguuuaca gc 22
 60
 <210> 72
 <211> 22
 <212> RNA
 <213> Homo sapiens
 65
 <400> 72
 cuuggcaccu agcaagcacu ca 22
 70
 <210> 73
 <211> 22

	<212> RNA	
	<213> Homo sapiens	
5	<400> 73 ugagaaccac gucugcucug ag	22
	<210> 74	
	<211> 22	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 74 uuauauaca accugauaag ug	22
15	<210> 75	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 75 ugagaccucu ggguucugag cu	22
	<210> 76	
	<211> 22	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 76 gcugacuccu aguccagggc uc	22
30	<210> 77	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 77 cuuucagucg gauguuugca gc	22
	<210> 78	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 78 uagcagcaca ucaugguuua ca	22
	<210> 79	
	<211> 23	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 79 agcuacauug ucugcugggu uuc	23
55	<210> 80	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	

	<400> 80	
	aggcaagaug cuggcauagc u	21
5	<210> 81	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 81	
	ucucacacag aaaucgcacc cgu	23
	<210> 82	
	<211> 22	
15	<212> RNA	
	<213> Homo sapiens	
	<400> 82	
	caucaucguc ucaaaugagu cu	22
20	<210> 83	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 83	
	ugaaggucua cugugugcca gg	22
	<210> 84	
30	<211> 17	
	<212> RNA	
	<213> Homo sapiens	
	<400> 84	
35	ucucgcuggg gccucca	17
	<210> 85	
	<211> 23	
	<212> RNA	
40	<213> Homo sapiens	
	<400> 85	
	uggaagacua gugauuuugu ugu	23
45	<210> 86	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
50	<400> 86	
	cagugcaaug augaaagggc au	22
	<210> 87	
	<211> 21	
55	<212> RNA	
	<213> Homo sapiens	
	<400> 87	
	cugaccuaug aaugacagc c	21

	<210> 88	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 88	
	uuguacaugg uaggcuuua uu	22
10		
	<210> 89	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 89	
	uucccuuugu cauccaugc cu	22
20		
	<210> 90	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 90	
	uucaaguaau ucaggauagg u	21
25		
	<210> 91	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 91	
	ucgaccggac cucgaccggc u	21
35		
	<210> 92	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
40		
	<400> 92	
	cuauacgacc ugcugccuuu cu	22
45		
	<210> 93	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 93	
	uuauaaagca augagacuga uu	22
50		
	<210> 94	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 94	
	ugugacuggu ugaccagagg gg	22
	<210> 95	
	<211> 23	

	<212> RNA	
	<213> Homo sapiens	
5	<400> 95 ucuuggagua ggucauuggg ugg	23
	<210> 96	
	<211> 22	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 96 uguaaacauc cuacacucag cu	22
15	<210> 97	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 97 aaaagcuggg uugagagggc ga	22
	<210> 98	
	<211> 22	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 98 caagcuugua ucuauaggua ug	22
30	<210> 99	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 99 ugcggggcua gggcuaacag ca	22
	<210> 100	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 100 accaucgacc guugauugua cc	22
	<210> 101	
	<211> 22	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 101 uggcaguguc uuagcugguu gu	22
55	<210> 102	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	

	<400> 102 accacugacc guugacugua cc	22
5	<210> 103 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 103 ugauauguuu gauauuuag gu	22
15	<210> 104 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 104 uaacagucua cagccauggu cg	22
25	<210> 105 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 105 aacauucaac euguegguga gu	22
35	<210> 106 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 106 uagcaccauc ugaaueggua ua	22
45	<210> 107 <211> 23 <212> RNA <213> Homo sapiens	
50	<400> 107 cagugcaaua guauugucua age	23
55	<210> 108 <211> 21 <212> RNA <213> Homo sapiens	
	<400> 108 uaguagaccg uauagcguac g	21
	<210> 109 <211> 21 <212> RNA <213> Homo sapiens	
	<400> 109 ucacagugaa ccggucucuu u	21

<210> 110
 <211> 17
 <212> RNA
 <213> Homo sapiens
 5
 <400> 110
 gggagaaggg ucggggc 17
 10
 <210> 111
 <211> 23
 <212> RNA
 <213> Homo sapiens
 15
 <400> 111
 aaugacacga ucacucccg uga 23
 20
 <210> 112
 <211> 21
 <212> RNA
 <213> Homo sapiens
 <400> 112
 acucuuuccc uguugcacua c 21
 25
 <210> 113
 <211> 22
 <212> RNA
 <213> Homo sapiens
 30
 <400> 113
 cagugcaaug uaaaaagggc au 22
 35
 <210> 114
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 114
 cuuucaguca gauguugcu gc 22
 40
 <210> 115
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 115
 uggugggccg cagaacaugu gc 22
 50
 <210> 116
 <211> 23
 <212> RNA
 <213> Homo sapiens
 55
 <400> 116
 caaagugcug uucgucagg uag 23
 <210> 117
 <211> 22

	<212> RNA	
	<213> Homo sapiens	
5	<400> 117 aaucguacag ggucauccac uu	22
	<210> 118	
	<211> 22	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 118 ucaggcucag ucccucccg au	22
15	<210> 119	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 119 uggcucaguu cagcaggaac ag	22
	<210> 120	
	<211> 22	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 120 ucugugagac caaagaacua cu	22
30	<210> 121	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 121 ucuggcuccg ugucuucacu ccc	23
	<210> 122	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 122 uucaccaccu ucuccacca gc	22
	<210> 123	
	<211> 23	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 123 uuuggcacua gcacauuuuu gcu	23
55	<210> 124	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	

	<400> 124 acucggcgug gcgucggucg ug	22
5	<210> 125 <211> 23 <212> RNA <213> Homo sapiens	
10	<400> 125 aggcagugua guuagcugau ugc	23
15	<210> 126 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 126 gccugcuggg guggaaccug gu	22
25	<210> 127 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 127 aaguucuguu auacacucag gc	22
35	<210> 128 <211> 23 <212> RNA <213> Homo sapiens	
40	<400> 128 ucaagagcaa uaacgaaaaa ugu	23
45	<210> 129 <211> 23 <212> RNA <213> Homo sapiens	
50	<400> 129 gauugagacu aguagggcua ggc	23
55	<210> 130 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 130 agggcuuagc ugcugugag ca	22
	<210> 131 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 131 aaugcacgg uauccaucug ua	22

<210> 132
 <211> 21
 <212> RNA
 <213> Homo sapiens
 5
 <400> 132
 ugucuugcag gccgucaugc a 21
 10
 <210> 133
 <211> 20
 <212> RNA
 <213> Homo sapiens
 15
 <400> 133
 guagaggaga uggcgcaggg 20
 20
 <210> 134
 <211> 23
 <212> RNA
 <213> Homo sapiens
 <400> 134
 agugccugag ggaguaagag ccc 23
 25
 <210> 135
 <211> 17
 <212> RNA
 <213> Homo sapiens
 30
 <400> 135
 gcggggcugg gcgcgcg 17
 35
 <210> 136
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 136
 uggugggcac agaaucugga cu 22
 40
 <210> 137
 <211> 23
 <212> RNA
 <213> Homo sapiens
 45
 <400> 137
 uauggcuuuu cauuccaug uga 23
 50
 <210> 138
 <211> 21
 <212> RNA
 <213> Homo sapiens
 55
 <400> 138
 uaccacaggg uagaaccacg g 21
 <210> 139
 <211> 24

	<212> RNA	
	<213> Homo sapiens	
5	<400> 139 aauccuugga accuaggugu gagu	24
	<210> 140	
	<211> 21	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 140 gcaguccaug ggcauauaca c	21
15	<210> 141	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 141 uuugugaccu gguccacuaa cc	22
	<210> 142	
	<211> 21	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 142 uacaguacug ugauaacuga a	21
30	<210> 143	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 143 auauauaca accugcuaag ug	22
	<210> 144	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 144 aaaguucuga gacacuccga cu	22
	<210> 145	
	<211> 23	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 145 caaagugcuu acagucagg uag	23
55	<210> 146	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	

	<400> 146 uaaagugcuu auagugcagg uag	23
5	<210> 147 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 147 cugcccuggc ccgagggacc ga	22
15	<210> 148 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 148 aacuggcccu caaagucccg cu	22
25	<210> 149 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 149 caaaaacugc aguuacuuuu gc	22
35	<210> 150 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 150 aucauacaag gacaauuucu uu	22
45	<210> 151 <211> 23 <212> RNA <213> Homo sapiens	
50	<400> 151 aucaacagac auuaauuggg cgc	23
55	<210> 152 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 152 aaggagcuca cagucuauug ag	22
	<210> 153 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 153 gucauacacg gcucuccucu cu	22

	<210> 154	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 154	
	ugaggcucug uuagccuugg cuc	23
10		
	<210> 155	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 155	
	cgucccgggg cugcgcgagg ca	22
20		
	<210> 156	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 156	
	ugccuacuga gcugaaacac ag	22
25		
	<210> 157	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 157	
	aacauucauu guugucggug ggu	23
35		
	<210> 158	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
40		
	<400> 158	
	cacauuacac ggucgaccuc u	21
45		
	<210> 159	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 159	
	ccgcacugug gguacuugcu gc	22
50		
	<210> 160	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 160	
	acaggugagg uucuugggag cc	22
	<210> 161	
	<211> 22	

<212> RNA
 <213> Homo sapiens

5 <400> 161
 ucucugggcc ugugucuag gc 22

10 <210> 162
 <211> 17
 <212> RNA
 <213> Homo sapiens

15 <400> 162
 gugggggaga ggcuguc 17

20 <210> 163
 <211> 23
 <212> RNA
 <213> Homo sapiens

25 <400> 163
 ugugcaaauc caugcaaac uga 23

30 <210> 164
 <211> 23
 <212> RNA
 <213> Homo sapiens

35 <400> 164
 cagugcaaug auauugucaa age 23

40 <210> 165
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 165
 agaggcuggc cgugaau uc 22

50 <210> 166
 <211> 23
 <212> RNA
 <213> Homo sapiens

55 <400> 166
 uagcaccau ugaauucagu guu 23

60 <210> 167
 <211> 22
 <212> RNA
 <213> Homo sapiens

65 <400> 167
 aaggguucc ucucugcagg ac 22

70 <210> 168
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 168	
	caggucgucu ugcagggcuu cu	22
5	<210> 169	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 169	
	uagugcaaua uugcuuauag ggu	23
	<210> 170	
	<211> 21	
	<212> RNA	
15	<213> Homo sapiens	
	<400> 170	
	uaaagugcug acagugcaga u	21
20	<210> 171	
	<211> 19	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 171	
	accgugcaaa gguagcaua	19
	<210> 172	
	<211> 22	
30	<212> RNA	
	<213> Homo sapiens	
	<400> 172	
35	ugcuaugcca acauauugcc au	22
	<210> 173	
	<211> 22	
	<212> RNA	
40	<213> Homo sapiens	
	<400> 173	
	cuuauagau uguauuguaa uu	22
45	<210> 174	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
50	<400> 174	
	aucaugaugg gcuccucggu gu	22
	<210> 175	
	<211> 18	
	<212> RNA	
55	<213> Homo sapiens	
	<400> 175	
	ggugggcuuc ccggaggg	18

	<210> 176	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 176	
	ugagaugaag cacuguagcu c	21
10		
	<210> 177	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 177	
	ugugcaaauc uaugcaaaac uga	23
20		
	<210> 178	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 178	
	caugccuuga guguaggacc gu	22
25		
	<210> 179	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 179	
	aucaaggau uaaaacuuug cc	22
35		
	<210> 180	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
40		
	<400> 180	
	gguggcccg ccgugccuga gg	22
45		
	<210> 181	
	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
	<400> 181	
	uugcagcugc cugggaguga cuuc	24
50		
	<210> 182	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 182	
	uauugggau gguaaaccgc uu	22
	<210> 183	
	<211> 22	

	<212> RNA	
	<213> Homo sapiens	
5	<400> 183 auaaagcuag auaaccgaaa gu	22
	<210> 184	
	<211> 22	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 184 acugcaguga aggcacuugu ag	22
15	<210> 185	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 185 aacauagagg aaauccacg u	21
	<210> 186	
	<211> 22	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 186 cagcagcaau ucauguuuug aa	22
30	<210> 187	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 187 uacccauugc auaucggagu ug	22
	<210> 188	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 188 uugcauaguc acaaaaguga uc	22
	<210> 189	
	<211> 23	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 189 ugaggaugga uagcaaggaa gcc	23
55	<210> 190	
	<211> 24	
	<212> RNA	
	<213> Homo sapiens	

<400> 190
 cccggacagg cguucgugcg acgu 24

 5 <210> 191
 <211> 18
 <212> RNA
 <213> Homo sapiens

 10 <400> 191
 gggcucacau caccuccau 18

 15 <210> 192
 <211> 22
 <212> RNA
 <213> Homo sapiens

 20 <400> 192
 uaugugccuu uggacuacau cg 22

 25 <210> 193
 <211> 22
 <212> RNA
 <213> Homo sapiens

 30 <400> 193
 aaacauucgc ggugcacuuc uu 22

 35 <210> 194
 <211> 20
 <212> RNA
 <213> Homo sapiens

 40 <400> 194
 uaaagagccc uguggagaca 20

 45 <210> 195
 <211> 23
 <212> RNA
 <213> Homo sapiens

 50 <400> 195
 gcaaagcaca cggccugcag aga 23

 55 <210> 196
 <211> 21
 <212> RNA
 <213> Homo sapiens

 <400> 196
 aaauaauacau gguugaucuu u 21

 <210> 197
 <211> 22
 <212> RNA
 <213> Homo sapiens

 <400> 197
 ugagaccagg acuggaugca cc 22

	<210> 198	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 198	
	aaaaguacuu gcggauuuug cu	22
10		
	<210> 199	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 199	
	ucaccagccc uguguuccu ag	22
20		
	<210> 200	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 200	
	cucccacaug caggguuugc a	21
25		
	<210> 201	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 201	
	agagcuuagc ugauugguga ac	22
35		
	<210> 202	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
40		
	<400> 202	
	agagguugcc cuuggugaau uc	22
45		
	<210> 203	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 203	
	aaucacuacag ggacauccag uu	22
50		
	<210> 204	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 204	
	agguugggau cgguugcaau gcu	23
	<210> 205	
	<211> 22	

<212> RNA
 <213> Homo sapiens

5 <400> 205
 auguagggcu aaaagccaug gg 22

10 <210> 206
 <211> 21
 <212> RNA
 <213> Homo sapiens

15 <400> 206
 uuugucuuga ucuaaccaug u 21

20 <210> 207
 <211> 23
 <212> RNA
 <213> Homo sapiens

25 <400> 207
 uagccccag gcuucacuug gcg 23

30 <210> 208
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 208
 agggacggga cgcgugcag ug 22

40 <210> 209
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 209
 auauacaggg ggagacucuu au 22

50 <210> 210
 <211> 21
 <212> RNA
 <213> Homo sapiens

55 <400> 210
 accacugcac uccagccuga g 21

60 <210> 211
 <211> 21
 <212> RNA
 <213> Homo sapiens

65 <400> 211
 auccccagau acauggaca a 21

70 <210> 212
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 212 gggguuccug gggaugggau uu	22
5	<210> 213 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 213 cugggagagg guuguuuacu cc	22
15	<210> 214 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 214 aacacaccug guuaaccucu uu	22
25	<210> 215 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 215 cuccuauaug augccuuucu uc	22
35	<210> 216 <211> 24 <212> RNA <213> Homo sapiens	
40	<400> 216 caaagugaug aguaauacug gcug	24
45	<210> 217 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 217 cuccugacuc cagguccugu gu	22
55	<210> 218 <211> 23 <212> RNA <213> Homo sapiens	
	<400> 218 gaggcugaug ugaguagacc acu	23
	<210> 219 <211> 23 <212> RNA <213> Homo sapiens	
	<400> 219 ugaguggggc uccgggacg gcg	23

<210> 220
 <211> 21
 <212> RNA
 <213> Homo sapiens
 5
 <400> 220
 guuucaccau guuggucagg c 21
 10
 <210> 221
 <211> 21
 <212> RNA
 <213> Homo sapiens
 15
 <400> 221
 aaauuuauac agucaaccuc u 21
 20
 <210> 222
 <211> 21
 <212> RNA
 <213> Homo sapiens
 <400> 222
 cuauacaauc uacugucuuu c 21
 25
 <210> 223
 <211> 22
 <212> RNA
 <213> Homo sapiens
 30
 <400> 223
 uagcagcaca uaaugguuug ug 22
 35
 <210> 224
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 224
 uggagagaaa ggcaguuccu ga 22
 40
 <210> 225
 <211> 21
 <212> RNA
 <213> Homo sapiens
 45
 <400> 225
 aggcggagac uugggcaauu g 21
 50
 <210> 226
 <211> 23
 <212> RNA
 <213> Homo sapiens
 55
 <400> 226
 uacaacaaau cacugaugcu gga 23
 <210> 227
 <211> 19

<212> RNA
 <213> Homo sapiens

5 <400> 227
 acuggccugg gacuaccgg 19

10 <210> 228
 <211> 23
 <212> RNA
 <213> Homo sapiens

15 <400> 228
 ugagcgccuc gacgacagag ccg 23

20 <210> 229
 <211> 22
 <212> RNA
 <213> Homo sapiens

25 <400> 229
 cuuagcaggu uguauuauca uu 22

30 <210> 230
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 230
 auggccagag cucacacaga gg 22

40 <210> 231
 <211> 20
 <212> RNA
 <213> Homo sapiens

45 <400> 231
 ggcuccuugg ucuaggggua 20

50 <210> 232
 <211> 17
 <212> RNA
 <213> Homo sapiens

55 <400> 232
 cuccgggacg gcugggc 17

60 <210> 233
 <211> 22
 <212> RNA
 <213> Homo sapiens

65 <400> 233
 gcuaaggaag uccugugcuc ag 22

70 <210> 234
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 234 ggagaaaua uccuuggugu gu	22
5	<210> 235 <211> 20 <212> RNA <213> Homo sapiens	
10	<400> 235 aaaagugauu gcaguguuug	20
15	<210> 236 <211> 21 <212> RNA <213> Homo sapiens	
20	<400> 236 ccaguccugu gccugccgcc u	21
25	<210> 237 <211> 21 <212> RNA <213> Homo sapiens	
30	<400> 237 aucauagagg aaaauccacg u	21
35	<210> 238 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 238 gaaguuguuc gugguggauu cg	22
45	<210> 239 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 239 cagcccggau cccagcccac uu	22
55	<210> 240 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 240 ugaaacauac acgggaaacc uc	22
	<210> 241 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 241 aaacaaacau ggugcacuuc uu	22

	<210> 242	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 242	
	ucagcaaaca uuauuugugu gc	22
10		
	<210> 243	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 243	
	caagcucgug ucuguggguc cg	22
20		
	<210> 244	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 244	
	uaggacacau ggucuacuuc u	21
25		
	<210> 245	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 245	
	cucacugaac aaugaaugca a	21
35		
	<210> 246	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
40		
	<400> 246	
	accuuggcuc uagacugcuu acu	23
45		
	<210> 247	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 247	
	caccuugcgc uacucagguc ug	22
50		
	<210> 248	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 248	
	uccgucucag uuacuuuaua gc	22
	<210> 249	
	<211> 22	

<212> RNA
 <213> Homo sapiens

5 <400> 249
 gcaugugaug aagcaauca gu 22

10 <210> 250
 <211> 23
 <212> RNA
 <213> Homo sapiens

15 <400> 250
 ucccccaggu gugauucuga uuu 23

20 <210> 251
 <211> 17
 <212> RNA
 <213> Homo sapiens

25 <400> 251
 ggcgggugcg ggggugg 17

30 <210> 252
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 252
 ccucccacac ccaaggcuug ca 22

40 <210> 253
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 253
 cacgcucaug cacacacca ca 22

50 <210> 254
 <211> 23
 <212> RNA
 <213> Homo sapiens

55 <400> 254
 agcagcauug uacagggcua uca 23

60 <210> 255
 <211> 22
 <212> RNA
 <213> Homo sapiens

65 <400> 255
 cugaagcuca gagggcucug au 22

70 <210> 256
 <211> 23
 <212> RNA
 <213> Homo sapiens

	<400> 256 uaaggugcau cuagugcaga uag	23
5	<210> 257 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 257 ccuauucuug auuacuuguu uc	22
15	<210> 258 <211> 21 <212> RNA <213> Homo sapiens	
20	<400> 258 agggccccc cucauccug u	21
25	<210> 259 <211> 21 <212> RNA <213> Homo sapiens	
30	<400> 259 agccgcgggg aucgccgagg g	21
35	<210> 260 <211> 21 <212> RNA <213> Homo sapiens	
40	<400> 260 aaucauucac ggacaacacu u	21
45	<210> 261 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 261 uacgcgcaga ccacaggaug uc	22
55	<210> 262 <211> 21 <212> RNA <213> Homo sapiens	
	<400> 262 cuggauggcu ccuccauguc u	21
	<210> 263 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 263 aguugccuuu uuguuccau gc	22

<210> 264
 <211> 18
 <212> RNA
 <213> Homo sapiens
 5
 <400> 264
 gggugcgggc cggcgggg 18
 10
 <210> 265
 <211> 22
 <212> RNA
 <213> Homo sapiens
 15
 <400> 265
 acccuaucaa uauugucucu gc 22
 20
 <210> 266
 <211> 23
 <212> RNA
 <213> Homo sapiens
 25
 <400> 266
 ccggucccag gagaaccugc aga 23
 30
 <210> 267
 <211> 22
 <212> RNA
 <213> Homo sapiens
 35
 <400> 267
 ugaguauuac auggccauc uc 22
 40
 <210> 268
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 268
 ccaaaacugc aguuacuuuu gc 22
 50
 <210> 269
 <211> 27
 <212> RNA
 <213> Homo sapiens
 55
 <400> 269
 accuucuugu auaagcacug ugcuaaa 27
 60
 <210> 270
 <211> 24
 <212> RNA
 <213> Homo sapiens
 65
 <400> 270
 agccuggaag cuggagccug cagu 24
 70
 <210> 271
 <211> 22

<212> RNA
 <213> Homo sapiens

5 <400> 271
 uuagggcccu ggcuccaucu cc 22

10 <210> 272
 <211> 23
 <212> RNA
 <213> Homo sapiens

15 <400> 272
 acuccauuug uuuugaugau gga 23

20 <210> 273
 <211> 23
 <212> RNA
 <213> Homo sapiens

25 <400> 273
 cccaguguuc agacuaccug uuc 23

30 <210> 274
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 274
 gagggauuggg uggaggcucu cc 22

40 <210> 275
 <211> 21
 <212> RNA
 <213> Homo sapiens

45 <400> 275
 ugcacggcac uggggacacg u 21

50 <210> 276
 <211> 20
 <212> RNA
 <213> Homo sapiens

55 <400> 276
 acugcccag gugcugcugg 20

60 <210> 277
 <211> 21
 <212> RNA
 <213> Homo sapiens

65 <400> 277
 gaacggcuuc auacaggagu u 21

70 <210> 278
 <211> 21
 <212> RNA
 <213> Homo sapiens

	<400> 278 aggggugcua ucugugauug a	21
5	<210> 279 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 279 uaaugccccc aaaaauccuu au	22
15	<210> 280 <211> 15 <212> RNA <213> Homo sapiens	
20	<400> 280 aggagauccu ggguu	15
25	<210> 281 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 281 aaugcaccug ggcaaggauu ca	22
35	<210> 282 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 282 cgucaacacu ugcugguuuc cu	22
45	<210> 283 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 283 ugucuacuc ccucaggcac au	22
55	<210> 284 <211> 21 <212> RNA <213> Homo sapiens	
	<400> 284 gccccgggca gugugaucu c	21
	<210> 285 <211> 24 <212> RNA <213> Homo sapiens	
	<400> 285 aaagacauag gauagaguca ccuc	24

<210> 286
 <211> 22
 <212> RNA
 <213> Homo sapiens
 5
 <400> 286
 auaauacaug guuaaccucu uu 22
 <210> 287
 <211> 23
 <212> RNA
 <213> Homo sapiens
 10
 <400> 287
 uauucauuua ucccagccu aca 23
 <210> 288
 <211> 23
 <212> RNA
 <213> Homo sapiens
 20
 <400> 288
 aggaagcccu ggaggggug gag 23
 <210> 289
 <211> 20
 <212> RNA
 <213> Homo sapiens
 25
 <400> 289
 cggcucuggg ucugugggga 20
 <210> 290
 <211> 22
 <212> RNA
 <213> Homo sapiens
 35
 <400> 290
 cuauacggcc uccuagcuu cc 22
 <210> 291
 <211> 18
 <212> RNA
 <213> Homo sapiens
 40
 <400> 291
 cgggcguggu gguggggg 18
 <210> 292
 <211> 19
 <212> RNA
 <213> Homo sapiens
 50
 <400> 292
 ggagauggag guugcagug 19
 <210> 293
 <211> 21

<212> RNA
 <213> Homo sapiens

5 <400> 293
 ugcaggacca agaugagccc u 21

10 <210> 294
 <211> 24
 <212> RNA
 <213> Homo sapiens

15 <400> 294
 uggcccugac ugaagaccag cagu 24

20 <210> 295
 <211> 22
 <212> RNA
 <213> Homo sapiens

25 <400> 295
 uaacacuguc ugguaaagau gg 22

30 <210> 296
 <211> 21
 <212> RNA
 <213> Homo sapiens

35 <400> 296
 cuccguuugc cuguuucgu g 21

40 <210> 297
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 297
 cuggcccucu cugcccuucc gu 22

50 <210> 298
 <211> 21
 <212> RNA
 <213> Homo sapiens

55 <400> 298
 caaacguga ggcgcugcua u 21

60 <210> 299
 <211> 20
 <212> RNA
 <213> Homo sapiens

65 <400> 299
 ggauccgagu cacggcacca 20

70 <210> 300
 <211> 17
 <212> RNA
 <213> Homo sapiens

	<400> 300	
	gagacugggg uggggcc	17
5	<210> 301	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 301	
	uuagugcaua gucuuugguc u	21
	<210> 302	
	<211> 22	
15	<212> RNA	
	<213> Homo sapiens	
	<400> 302	
	uuaauuuuuu guuucgguca cu	22
20	<210> 303	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 303	
	uaauccuugc uaccugggug aga	23
	<210> 304	
	<211> 22	
30	<212> RNA	
	<213> Homo sapiens	
	<400> 304	
	aaaaguaauu gugguuuugg cc	22
35	<210> 305	
	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
40	<400> 305	
	cugaagugau guguaacuga ucag	24
	<210> 306	
45	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 306	
50	auucuaauuu cuccacgucu uu	22
	<210> 307	
	<211> 22	
	<212> RNA	
55	<213> Homo sapiens	
	<400> 307	
	gacuauagaa cuuucccccu ca	22

<210> 308
 <211> 21
 <212> RNA
 <213> Homo sapiens
 5
 <400> 308
 aauggcgcca cuagguugu g 21
 10
 <210> 309
 <211> 20
 <212> RNA
 <213> Homo sapiens
 15
 <400> 309
 accaggaggc ugaggcccu 20
 20
 <210> 310
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 310
 acuccagccc cacagccuca gc 22
 25
 <210> 311
 <211> 23
 <212> RNA
 <213> Homo sapiens
 30
 <400> 311
 ccaguuaccg cuuccgcuac cgc 23
 35
 <210> 312
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 312
 auccgcguc ugacucucug cc 22
 40
 <210> 313
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 313
 uuuccggcuc gcgugggugu gu 22
 50
 <210> 314
 <211> 22
 <212> RNA
 <213> Homo sapiens
 55
 <400> 314
 auauacaggg ggagacucuc au 22
 <210> 315
 <211> 22

<212> RNA
 <213> Homo sapiens

5 <400> 315
 accguggcuu ucgauuguua cu 22

<210> 316
 <211> 22
 <212> RNA
 10 <213> Homo sapiens

<400> 316
 ccaauauuac ugugcugcuu ua 22

15 <210> 317
 <211> 23
 <212> RNA
 <213> Homo sapiens

20 <400> 317
 caaagugcuc auagucagg uag 23

<210> 318
 <211> 21
 25 <212> RNA
 <213> Homo sapiens

<400> 318
 ccucccaugc caagaacucc c 21

30 <210> 319
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 319
 ugguuuaccg ucccacauac au 22

40 <210> 320
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 320
 cugggaggug gauguuuacu uc 22

<210> 321
 <211> 22
 <212> RNA
 50 <213> Homo sapiens

<400> 321
 cugggagaag gcuguuuacu cu 22

55 <210> 322
 <211> 20
 <212> RNA
 <213> Homo sapiens

	<400> 322	
	uuggccaugg ggcugcgcg	20
5	<210> 323	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 323	
	ucucucggcu ccucgcggu c	21
	<210> 324	
	<211> 22	
15	<212> RNA	
	<213> Homo sapiens	
	<400> 324	
	ucaccugca uccgcaccc ag	22
20	<210> 325	
	<211> 19	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 325	
	gacuggacaa gcugaggaa	19
	<210> 326	
30	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
	<400> 326	
35	gaaaaugaug aguagugacu gaug	24
	<210> 327	
	<211> 22	
	<212> RNA	
40	<213> Homo sapiens	
	<400> 327	
	uaguggauga ugcacucugu gc	22
45	<210> 328	
	<211> 17	
	<212> RNA	
	<213> Homo sapiens	
50	<400> 328	
	acccacucc ugguacc	17
	<210> 329	
	<211> 22	
55	<212> RNA	
	<213> Homo sapiens	
	<400> 329	
	cacccccugu uuccugccc ac	22

	<210> 330	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 330	
	acacaugggu ggcuguggcc u	21
10		
	<210> 331	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 331	
	ugugacagau ugauaacuga aa	22
20		
	<210> 332	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 332	
	cagggaaaug ggaagaacua ga	22
25		
	<210> 333	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 333	
	cgaaaacagc aaauaccuuu gc	22
35		
	<210> 334	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
40		
	<400> 334	
	ugagugugug ugugugagug ugu	23
45		
	<210> 335	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 335	
	ucuaguaaga guggcagucg a	21
50		
	<210> 336	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 336	
	uagucugcu gaccaucacc uu	22
	<210> 337	
	<211> 23	

<212> RNA
 <213> Homo sapiens

5 <400> 337
 cugggaucuc cggggucuug guu 23

<210> 338
 <211> 21
 <212> RNA
 10 <213> Homo sapiens

<400> 338
 cugacuguug ccguccucca g 21

15 <210> 339
 <211> 22
 <212> RNA
 <213> Homo sapiens

20 <400> 339
 cuauacaacc uacugccuuc cc 22

<210> 340
 <211> 26
 25 <212> RNA
 <213> Homo sapiens

<400> 340
 30 aaguaguugg uuuguauagag augguu 26

<210> 341
 <211> 23
 <212> RNA
 <213> Homo sapiens
 35

<400> 341
 aggaugagca aagaaaguag auu 23

40 <210> 342
 <211> 22
 <212> RNA
 <213> Homo sapiens

<400> 342
 45 uugcuugaac ccaggaagug ga 22

<210> 343
 <211> 23
 <212> RNA
 50 <213> Homo sapiens

<400> 343
 uggaguccag gaaucugcau uuu 23

55 <210> 344
 <211> 21
 <212> RNA
 <213> Homo sapiens

	<400> 344 ucagugcaug acagaacuug g	21
5	<210> 345 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 345 uguaacagca acuccaugug ga	22
15	<210> 346 <211> 21 <212> RNA <213> Homo sapiens	
20	<400> 346 uagcagcaca gaaauauugg c	21
25	<210> 347 <211> 23 <212> RNA <213> Homo sapiens	
30	<400> 347 uaauacugcc ggguaaugau gga	23
35	<210> 348 <211> 21 <212> RNA <213> Homo sapiens	
40	<400> 348 uaacagucuc cagucacggc c	21
45	<210> 349 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 349 cucaguagcc aguguagauc cu	22
55	<210> 350 <211> 23 <212> RNA <213> Homo sapiens	
	<400> 350 ucagcaccag gauauuguug gag	23
	<210> 351 <211> 20 <212> RNA <213> Homo sapiens	
	<400> 351 auauggguuu acuaguuggu	20

	<210> 352	
	<211> 22	
	<212> RNA	
5	<213> Homo sapiens	
	<400> 352	
	ugagggacag augccagaag ca	22
10	<210> 353	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
15	<400> 353	
	uagugaguua gagaugcaga gcc	23
20	<210> 354	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
	<400> 354	
	cgcauccccu agggc auugg ugu	23
25	<210> 355	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
30	<400> 355	
	gugcauugua guugcauugc a	21
35	<210> 356	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
40	<400> 356	
	cucgugggcu cuggccacgg cc	22
45	<210> 357	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 357	
	agaucgaccg uguuauauuc gc	22
50	<210> 358	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55	<400> 358	
	aucgggaug ucguguccgc cc	22
	<210> 359	
	<211> 17	

<212> RNA
 <213> Homo sapiens

5 <400> 359
 gaagauggac guacuuu 17

10 <210> 360
 <211> 22
 <212> RNA
 <213> Homo sapiens

15 <400> 360
 uggcggcggu aguuaugggc uu 22

20 <210> 361
 <211> 22
 <212> RNA
 <213> Homo sapiens

25 <400> 361
 uuucuauuuc ucaguggggc uc 22

30 <210> 362
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 362
 aggacuggac ucccggcagc cc 22

40 <210> 363
 <211> 18
 <212> RNA
 <213> Homo sapiens

45 <400> 363
 cggugagcgc ugcuggc 18

50 <210> 364
 <211> 22
 <212> RNA
 <213> Homo sapiens

55 <400> 364
 aggaccuucc cugaaccaag ga 22

60 <210> 365
 <211> 22
 <212> RNA
 <213> Homo sapiens

65 <400> 365
 aggcggggcg ccgcgggacc gc 22

70 <210> 366
 <211> 21
 <212> RNA
 <213> Homo sapiens

	<400> 366	
	aaggcagggc ccccgcuccc c	21
5	<210> 367	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 367	
	caagcucgcu ucuauggguc ug	22
	<210> 368	
	<211> 21	
15	<212> RNA	
	<213> Homo sapiens	
	<400> 368	
	agaggauacc cuuuguaugu u	21
20	<210> 369	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 369	
	ugagccccug ugccgcccc ag	22
	<210> 370	
30	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 370	
35	uccuucugcu ccguccccca g	21
	<210> 371	
	<211> 22	
	<212> RNA	
40	<213> Homo sapiens	
	<400> 371	
	agaaggaaau ugaaaucauu ua	22
45	<210> 372	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
50	<400> 372	
	agugaaugau gguucugac c	21
	<210> 373	
	<211> 19	
55	<212> RNA	
	<213> Homo sapiens	
	<400> 373	
	aucccaccac ugccaccau	19

<210> 374
 <211> 25
 <212> RNA
 <213> Homo sapiens
 5
 <400> 374
 gaacccauga gguugaggcu gcagu 25
 10
 <210> 375
 <211> 19
 <212> RNA
 <213> Homo sapiens
 15
 <400> 375
 uggauuuuug gaucaggga 19
 20
 <210> 376
 <211> 18
 <212> RNA
 <213> Homo sapiens
 <400> 376
 acguuggcuc ugguggug 18
 25
 <210> 377
 <211> 18
 <212> RNA
 <213> Homo sapiens
 30
 <400> 377
 cagggaggug aaugugau 18
 35
 <210> 378
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 378
 cuccuggggc ccgcacucuc gc 22
 40
 <210> 379
 <211> 23
 <212> RNA
 <213> Homo sapiens
 45
 <400> 379
 agcugguguu gugaaucagg ccg 23
 50
 <210> 380
 <211> 22
 <212> RNA
 <213> Homo sapiens
 55
 <400> 380
 cagugguuuu acccauggu ag 22
 <210> 381
 <211> 22

<212> RNA
 <213> Homo sapiens

5 <400> 381
 ugcccugugg acucaguucu gg 22

<210> 382
 <211> 22
 <212> RNA
 10 <213> Homo sapiens

<400> 382
 gcccaaaggu gaauuuuuug gg 22

15 <210> 383
 <211> 21
 <212> RNA
 <213> Homo sapiens

20 <400> 383
 cggcggggac ggcgauuggu c 21

25 <210> 384
 <211> 20
 <212> RNA
 <213> Homo sapiens

30 <400> 384
 ccccagggcg acgcggcggg 20

35 <210> 385
 <211> 22
 <212> RNA
 <213> Homo sapiens

<400> 385
 accuugccuu gcugcccggg cc 22

40 <210> 386
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 386
 aacuggccua caaaguucca gu 22

50 <210> 387
 <211> 23
 <212> RNA
 <213> Homo sapiens

<400> 387
 aguuuugcag guuugcaucc agc 23

55 <210> 388
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 388 auaagacgaa caaaagguuu gu	22
5	<210> 389 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 389 uuggggaaac ggccgcugag ug	22
15	<210> 390 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 390 agaguugagu cuggacgucc cg	22
25	<210> 391 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 391 ccuguucucc auuacuuggc uc	22
35	<210> 392 <211> 23 <212> RNA <213> Homo sapiens	
40	<400> 392 agaauugcgu uuggacaauc agu	23
45	<210> 393 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 393 acugauuucu uuugguguuc ag	22
55	<210> 394 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 394 caucuggcau ccgucacaca ga	22
	<210> 395 <211> 21 <212> RNA <213> Homo sapiens	
	<400> 395 gcugcaccgg agacugggua a	21

<210> 396
 <211> 21
 <212> RNA
 <213> Homo sapiens
 5
 <400> 396
 uaccagucu ccgugcagc c 21
 <210> 397
 <211> 21
 <212> RNA
 <213> Homo sapiens
 10
 <400> 397
 accugaauua ccaaagcuu u 21
 <210> 398
 <211> 21
 <212> RNA
 <213> Homo sapiens
 20
 <400> 398
 ccaggcucug caguggaac u 21
 <210> 399
 <211> 22
 <212> RNA
 <213> Homo sapiens
 25
 <400> 399
 cugcccuagu cuagcugaag cu 22
 <210> 400
 <211> 22
 <212> RNA
 <213> Homo sapiens
 30
 <400> 400
 uggggcggag cuuccggagg cc 22
 <210> 401
 <211> 23
 <212> RNA
 <213> Homo sapiens
 35
 <400> 401
 agguuguccg uggugaguuc gca 23
 <210> 402
 <211> 23
 <212> RNA
 <213> Homo sapiens
 40
 <400> 402
 ucccuguccu ccaggagcuc acg 23
 <210> 403
 <211> 22

<212> RNA
 <213> Homo sapiens

5 <400> 403
 caaucagcaa guauacugcc cu 22

10 <210> 404
 <211> 22
 <212> RNA
 <213> Homo sapiens

15 <400> 404
 caaucacuaa cuccacugcc au 22

20 <210> 405
 <211> 22
 <212> RNA
 <213> Homo sapiens

25 <400> 405
 aaucacuaac cacacggcca gg 22

30 <210> 406
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 406
 uuuaagaaaa caccauggag au 22

40 <210> 407
 <211> 23
 <212> RNA
 <213> Homo sapiens

45 <400> 407
 agggacuuuu gggggcagau gug 23

50 <210> 408
 <211> 20
 <212> RNA
 <213> Homo sapiens

55 <400> 408
 ccguguuucc cccacgcuuu 20

60 <210> 409
 <211> 23
 <212> RNA
 <213> Homo sapiens

65 <400> 409
 aguggaugau ggagacucgg uac 23

70 <210> 410
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 410	
	guagauucuc cuucuaugag ua	22
5	<210> 411	
	<211> 20	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 411	
	acugggcuug gagucagaag	20
	<210> 412	
	<211> 22	
15	<212> RNA	
	<213> Homo sapiens	
	<400> 412	
	uguccucuag ggccugcagu cu	22
20	<210> 413	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 413	
	ggaggaaccu ugagcuucg gc	22
	<210> 414	
30	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 414	
35	uuucagauaa caguauuaca u	21
	<210> 415	
	<211> 21	
	<212> RNA	
40	<213> Homo sapiens	
	<400> 415	
	ugugcagcag gccaacccag a	21
	<210> 416	
45	<211> 20	
	<212> RNA	
	<213> Homo sapiens	
	<400> 416	
50	ggcggcggcg gaggcggggg	20
	<210> 417	
	<211> 20	
55	<212> RNA	
	<213> Homo sapiens	
	<400> 417	
	uguuccucug ucucccagac	20

<210> 418
 <211> 18
 <212> RNA
 <213> Homo sapiens
 5
 <400> 418
 cucgaguugg aagaggcg 18
 10
 <210> 419
 <211> 22
 <212> RNA
 <213> Homo sapiens
 15
 <400> 419
 uggggauuug gagaaguggu ga 22
 20
 <210> 420
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 420
 auggcaucgu cccuggugg cu 22
 25
 <210> 421
 <211> 23
 <212> RNA
 <213> Homo sapiens
 30
 <400> 421
 agggaaaaaa aaaaggauuu guc 23
 35
 <210> 422
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 422
 uccaggcagg agccggacug ga 22
 40
 <210> 423
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 423
 uaggggcagc agaggaccug gg 22
 50
 <210> 424
 <211> 26
 <212> RNA
 <213> Homo sapiens
 55
 <400> 424
 cacaggacug acuccucacc ccagug 26
 <210> 425
 <211> 22

<212> RNA
 <213> Homo sapiens

5 <400> 425
 cacacaagug gcccccaaca cu 22

10 <210> 426
 <211> 20
 <212> RNA
 <213> Homo sapiens

15 <400> 426
 cgccccuccu gccccacag 20

20 <210> 427
 <211> 24
 <212> RNA
 <213> Homo sapiens

25 <400> 427
 ggugggaugg agagaaggua ugag 24

30 <210> 428
 <211> 21
 <212> RNA
 <213> Homo sapiens

35 <400> 428
 aguggaccga ggaaggaagg a 21

40 <210> 429
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 429
 aguggggaac ccuuccauga gg 22

50 <210> 430
 <211> 22
 <212> RNA
 <213> Homo sapiens

55 <400> 430
 aauccuuugu ccuuggguga ga 22

60 <210> 431
 <211> 22
 <212> RNA
 <213> Homo sapiens

65 <400> 431
 aauccacgcu gagcuuggca uc 22

70 <210> 432
 <211> 25
 <212> RNA
 <213> Homo sapiens

	<400> 432	
	aaaggauucu gcugucgguc ccacu	25
5	<210> 433	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 433	
	ucggggauca ucaugucacg aga	23
	<210> 434	
	<211> 21	
15	<212> RNA	
	<213> Homo sapiens	
	<400> 434	
	gcgaccuaa cuugguuca g	21
20	<210> 435	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 435	
	ucagcuacua ccucuauuag g	21
	<210> 436	
30	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
	<400> 436	
35	uagauaaaau auugguaccu g	21
	<210> 437	
	<211> 22	
	<212> RNA	
40	<213> Homo sapiens	
	<400> 437	
	ucaguuccag gccaccagg cu	22
	<210> 438	
45	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
	<400> 438	
50	ucagaacaaa ugccgguucc caga	24
	<210> 439	
	<211> 22	
	<212> RNA	
55	<213> Homo sapiens	
	<400> 439	
	acucaaaacc cuucagugac uu	22

	<210> 440	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 440	
	augcugacau auuuacuaga gg	22
10		
	<210> 441	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 441	
	uggguuuacg uugggagaac u	21
20		
	<210> 442	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 442	
	uucuuuugcc ucccagccua ca	22
25		
	<210> 443	
	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 443	
	acuggcuagg gaaaugauu ggau	24
35		
	<210> 444	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
40		
	<400> 444	
	gcagcagaga auaggacuac guc	23
45		
	<210> 445	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 445	
	uauguaauau gguccauc uu	22
50		
	<210> 446	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 446	
	uauguaacau gguccacuaa cu	22
	<210> 447	
	<211> 22	

	<212> RNA	
	<213> Homo sapiens	
5	<400> 447 cgaaucauua uuugcugcuc ua	22
	<210> 448	
	<211> 23	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 448 uuaaugcuua ucgugauagg ggu	23
15	<210> 449	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 449 aucacacaaa ggcaacuuuu gu	22
	<210> 450	
	<211> 21	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 450 gccccugggc cuauccuaga a	21
30	<210> 451	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 451 augaccuaug aaugacaga c	21
	<210> 452	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 452 aucauagagg aaaauccaug uu	22
	<210> 453	
	<211> 21	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 453 cggggcagcu caguacagga u	21
55	<210> 454	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	

	<400> 454 acugcugagc uagcacuucc cg	22
5	<210> 455 <211> 18 <212> RNA <213> Homo sapiens	
10	<400> 455 ucgaggagcu cacagucu	18
15	<210> 456 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 456 gugaacgggc gccaucccga gg	22
25	<210> 457 <211> 23 <212> RNA <213> Homo sapiens	
30	<400> 457 ccucagggcgu guagaacagg gcu	23
35	<210> 458 <211> 23 <212> RNA <213> Homo sapiens	
40	<400> 458 uaagugcuuc cauguuuuag uag	23
45	<210> 459 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 459 aaaaacugag acuacuuuug ca	22
55	<210> 460 <211> 23 <212> RNA <213> Homo sapiens	
	<400> 460 uaagugcuuc cauguuugag ugu	23
	<210> 461 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 461 uacgucaucg uugucaucgu ca	22

	<210> 462	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 462	
	cauuauuacu uuugguacgc g	21
10		
	<210> 463	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 463	
	uaauuuuauug uauaagcuag u	21
20		
	<210> 464	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 464	
	ucugggcaac aaagugagac cu	22
25		
	<210> 465	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 465	
	uaccucuguag aaccgaauuu gug	23
35		
	<210> 466	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
	<400> 466	
	cugauaagaa cagaggccca gau	23
40		
	<210> 467	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
45		
	<400> 467	
	agaauugugg cuggacauu gu	22
50		
	<210> 468	
	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 468	
	agcgcgggcu gagcgugcc aguc	24
	<210> 469	
	<211> 22	

<212> RNA
 <213> Homo sapiens

5 <400> 469
 uguuguacuu uuuuuuuugu uc 22

<210> 470
 <211> 22
 <212> RNA
 10 <213> Homo sapiens

<400> 470
 uaacugguug aacaacugaa cc 22

15 <210> 471
 <211> 23
 <212> RNA
 <213> Homo sapiens

20 <400> 471
 aaaagugcuu acagugcagg uag 23

<210> 472
 <211> 22
 25 <212> RNA
 <213> Homo sapiens

<400> 472
 aaucauacac gguugaccua uu 22

30 <210> 473
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 473
 uagcaccauu ugaaaucggu ua 22

40 <210> 474
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 474
 gcugcguug gauuucgucc cc 22

<210> 475
 <211> 22
 <212> RNA
 50 <213> Homo sapiens

<400> 475
 uaauacugcc ugguaaugau ga 22

55 <210> 476
 <211> 23
 <212> RNA
 <213> Homo sapiens

	<400> 476 ugucacucgg cucggccac uac	23
5	<210> 477 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 477 ugagaacuga auuccauggg uu	22
15	<210> 478 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 478 cugugcgugu gacagcggcu ga	22
25	<210> 479 <211> 21 <212> RNA <213> Homo sapiens	
30	<400> 479 cagcagcaca cugugguuug u	21
35	<210> 480 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 480 caacaaauca cagucugcca ua	22
45	<210> 481 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 481 ugaccuggga cucggacagc ug	22
55	<210> 482 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 482 uauguaacac gguccacuaa cc	22
	<210> 483 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 483 cugcgcaagc uacugccuug cu	22

	<210> 484	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 484	
	aacuagcucu guggauccug ac	22
10		
	<210> 485	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 485	
	uuauuguuug ccugggacug ag	22
20		
	<210> 486	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 486	
	cuguugccac uaaccucaac cu	22
25		
	<210> 487	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 487	
	uccaguacca cgugucaggg cca	23
35		
	<210> 488	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 488	
	uccugucuuu ccuuguugga gc	22
40		
	<210> 489	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
45		
	<400> 489	
	uuacaguugu ucaaccaguu acu	23
50		
	<210> 490	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 490	
	acgcccucc ccccuucuu ca	22
	<210> 491	
	<211> 22	

	<212> RNA	
	<213> Homo sapiens	
5	<400> 491 gcucugacuu uauugcacua cu	22
	<210> 492	
	<211> 22	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 492 uucugccucu guccaggucc uu	22
15	<210> 493	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 493 aaucagugaa ugccuugaac cu	22
	<210> 494	
	<211> 22	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 494 acaguagucu gcacauuggu ua	22
30	<210> 495	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 495 aguucuucag uggcaagcuu ua	22
	<210> 496	
40	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
	<400> 496 ccuccguguu accuguccuc uag	23
	<210> 497	
	<211> 22	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 497 agaccuggu cugcacucua uc	22
55	<210> 498	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	

	<400> 498	
	uucuggaauu cugugugagg ga	22
5	<210> 499	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 499	
	caucccuugc augguggagg g	21
	<210> 500	
	<211> 22	
15	<212> RNA	
	<213> Homo sapiens	
	<400> 500	
	cuagguaugg ucccagggau cc	22
20	<210> 501	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 501	
	uagguuaucc guguugccuu cg	22
	<210> 502	
30	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
	<400> 502	
35	acugcccuaa gugcuccuuc ugg	23
	<210> 503	
	<211> 23	
	<212> RNA	
40	<213> Homo sapiens	
	<400> 503	
	uaagugcuuc cauguuuugg uga	23
45	<210> 504	
	<211> 20	
	<212> RNA	
	<213> Homo sapiens	
50	<400> 504	
	ccucugggcc cuuccuccag	20
	<210> 505	
	<211> 21	
55	<212> RNA	
	<213> Homo sapiens	
	<400> 505	
	cuguacaggc cacugccuug c	21

	<210> 506		
	<211> 22		
	<212> RNA		
	<213> Homo sapiens		
5			
	<400> 506		
	ugccuacuga gcugauauca gu	22	
10			
	<210> 507		
	<211> 25		
	<212> RNA		
	<213> Homo sapiens		
15			
	<400> 507		
	aaauaugaug aaacucacag cugag	25	
20			
	<210> 508		
	<211> 22		
	<212> RNA		
	<213> Homo sapiens		
	<400> 508		
	aagaugugga aaaauggaa uc	22	
25			
	<210> 509		
	<211> 23		
	<212> RNA		
	<213> Homo sapiens		
30			
	<400> 509		
	cucucaccac ugcccuccca cag	23	
35			
	<210> 510		
	<211> 22		
	<212> RNA		
	<213> Homo sapiens		
	<400> 510		
	uggguuccug gcaugcugau uu	22	
40			
	<210> 511		
	<211> 23		
	<212> RNA		
	<213> Homo sapiens		
45			
	<400> 511		
	acuuaaacgu ggauguacuu gcu	23	
50			
	<210> 512		
	<211> 22		
	<212> RNA		
	<213> Homo sapiens		
55			
	<400> 512		
	aaaucccuug uagauaaccc gg	22	
	<210> 513		
	<211> 22		

<212> RNA
 <213> Homo sapiens

5 <400> 513
 cuguacagcc uccuagcuuu cc 22

<210> 514
 <211> 22
 <212> RNA
 10 <213> Homo sapiens

<400> 514
 caaaauucgua ucuaggggaa ua 22

15 <210> 515
 <211> 22
 <212> RNA
 <213> Homo sapiens

20 <400> 515
 ugcuggauca gugguucgag uc 22

<210> 516
 <211> 23
 25 <212> RNA
 <213> Homo sapiens

<400> 516
 guccaguuuu cccaggauc ccu 23

30 <210> 517
 <211> 24
 <212> RNA
 <213> Homo sapiens

35 <400> 517
 gcugguuuca uauggugguu uaga 24

40 <210> 518
 <211> 23
 <212> RNA
 <213> Homo sapiens

45 <400> 518
 ucuggcaagu aaaaaacucu cau 23

<210> 519
 <211> 20
 <212> RNA
 50 <213> Homo sapiens

<400> 519
 gugcauugcu guugcauugc 20

55 <210> 520
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 520 uuuccuaccc uaccugaaga cu	22
5	<210> 521 <211> 21 <212> RNA <213> Homo sapiens	
10	<400> 521 acagggccgc agauggagac u	21
15	<210> 522 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 522 gaaaucaagc gugggugaga cc	22
25	<210> 523 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 523 uuugaggcua cagugagaug ug	22
35	<210> 524 <211> 23 <212> RNA <213> Homo sapiens	
40	<400> 524 cggcccgggc ugcugcuguu ccu	23
45	<210> 525 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 525 aaccaucgac cguugagugg ac	22
55	<210> 526 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 526 ugggucuuug cgggcgagau ga	22
	<210> 527 <211> 21 <212> RNA <213> Homo sapiens	
	<400> 527 ugauugucca aacgcaauuc u	21

	<210> 528	
	<211> 22	
	<212> RNA	
5	<213> Homo sapiens	
	<400> 528	
	uaggauuaca agugucggcc ac	22
10	<210> 529	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15	<400> 529	
	ugcccugccu guuuucuccu uu	22
20	<210> 530	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 530	
	cggggagaga acgcagugac gu	22
25	<210> 531	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
30	<400> 531	
	aaucugagaa ggcgcacaag gu	22
35	<210> 532	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
40	<400> 532	
	uugaagagga ggugcucugu agc	23
45	<210> 533	
	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
	<400> 533	
	gaugcgccgc ccacugcccc gcgc	24
50	<210> 534	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
55	<400> 534	
	uuaagacuug cagugauguu u	21
	<210> 535	
	<211> 23	

<212> RNA
 <213> Homo sapiens

5 <400> 535
 caaagacugc aaauacuuuu gcg 23

10 <210> 536
 <211> 25
 <212> RNA
 <213> Homo sapiens

15 <400> 536
 uggaacggg uuccggcaga cgcug 25

20 <210> 537
 <211> 22
 <212> RNA
 <213> Homo sapiens

25 <400> 537
 ccugugccc ggcccacuuc ug 22

30 <210> 538
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 538
 ccaauauugg cugugcugcu cc 22

40 <210> 539
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 539
 uuuaacaugg gguuaccugc ug 22

50 <210> 540
 <211> 22
 <212> RNA
 <213> Homo sapiens

55 <400> 540
 caaccugac gaucuccuca gc 22

60 <210> 541
 <211> 22
 <212> RNA
 <213> Homo sapiens

65 <400> 541
 cccaauacac ggucgaccuc uu 22

70 <210> 542
 <211> 22
 <212> RNA
 <213> Homo sapiens

<400> 542
 uuuuucauua uugcuccuga cc 22

 5 <210> 543
 <211> 20
 <212> RNA
 <213> Homo sapiens

 10 <400> 543
 acuguaaacg cuuucugaug 20

 15 <210> 544
 <211> 18
 <212> RNA
 <213> Homo sapiens

 20 <400> 544
 cuaagaaguu gacugaag 18

 25 <210> 545
 <211> 23
 <212> RNA
 <213> Homo sapiens

 30 <400> 545
 ugagcaccac acaggccggg cgc 23

 35 <210> 546
 <211> 21
 <212> RNA
 <213> Homo sapiens

 40 <400> 546
 auaggcacca aaaagcaaca a 21

 45 <210> 547
 <211> 22
 <212> RNA
 <213> Homo sapiens

 50 <400> 547
 uggcagugua uuguuagcug gu 22

 55 <210> 548
 <211> 22
 <212> RNA
 <213> Homo sapiens

 <400> 548
 gaagaacugu ugcauuugcc cu 22

 <210> 549
 <211> 19
 <212> RNA
 <213> Homo sapiens

 <400> 549
 cauccguccg ucuguccac 19

<210> 550
 <211> 22
 <212> RNA
 <213> Homo sapiens
 5
 <400> 550
 gggagccagg aaguauugau gu 22
 10
 <210> 551
 <211> 22
 <212> RNA
 <213> Homo sapiens
 15
 <400> 551
 caaacuggc aauuacuuuu gc 22
 20
 <210> 552
 <211> 21
 <212> RNA
 <213> Homo sapiens
 <400> 552
 uccucuucuc ccuccucca g 21
 25
 <210> 553
 <211> 23
 <212> RNA
 <213> Homo sapiens
 30
 <400> 553
 agcuucuuua cagugcugcc uug 23
 35
 <210> 554
 <211> 20
 <212> RNA
 <213> Homo sapiens
 <400> 554
 cgggcguggu gguggggug 20
 40
 <210> 555
 <211> 23
 <212> RNA
 <213> Homo sapiens
 45
 <400> 555
 cuggagauau ggaagagcug ugu 23
 50
 <210> 556
 <211> 22
 <212> RNA
 <213> Homo sapiens
 55
 <400> 556
 uggguggucu ggagauuugu gc 22
 <210> 557
 <211> 19

<212> RNA
 <213> Homo sapiens

5 <400> 557
 gaugaugcug cugaugcug 19

10 <210> 558
 <211> 22
 <212> RNA
 <213> Homo sapiens

15 <400> 558
 ucucccaacc cuuguaccag ug 22

20 <210> 559
 <211> 21
 <212> RNA
 <213> Homo sapiens

25 <400> 559
 ugauauguuu gauauugggu u 21

30 <210> 560
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 560
 cgggguuuug agggcgagau ga 22

40 <210> 561
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 561
 acugcauuau gagcacuuaa ag 22

50 <210> 562
 <211> 22
 <212> RNA
 <213> Homo sapiens

55 <400> 562
 uaaucucagc uggcaacugu ga 22

60 <210> 563
 <211> 23
 <212> RNA
 <213> Homo sapiens

65 <400> 563
 caggcaguga cuguucagac guc 23

70 <210> 564
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 564 agauguccag ccacaaucu cg	22
5	<210> 565 <211> 23 <212> RNA <213> Homo sapiens	
10	<400> 565 uguguacaca cgugccaggc gcu	23
15	<210> 566 <211> 20 <212> RNA <213> Homo sapiens	
20	<400> 566 aaaagcuggg uugagaggu	20
25	<210> 567 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 567 aggugguccg uggcgcuuc gc	22
35	<210> 568 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 568 caggcacggg agcucaggug ag	22
45	<210> 569 <211> 21 <212> RNA <213> Homo sapiens	
50	<400> 569 uaacgcauaa uauggacaug u	21
55	<210> 570 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 570 ucaggugugg aaacugaggc ag	22
	<210> 571 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 571 aagcaauacu guuaccugaa au	22

<210> 572
 <211> 21
 <212> RNA
 <213> Homo sapiens
 5
 <400> 572
 gaccgagagg gccucggcug u 21
 10
 <210> 573
 <211> 23
 <212> RNA
 <213> Homo sapiens
 15
 <400> 573
 ugggccaggg agcagcuggu ggg 23
 20
 <210> 574
 <211> 22
 <212> RNA
 <213> Homo sapiens
 25
 <400> 574
 accgaagacu gugcguaau cu 22
 30
 <210> 575
 <211> 22
 <212> RNA
 <213> Homo sapiens
 35
 <400> 575
 acaacaguga cuugcucucc aa 22
 40
 <210> 576
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 576
 ugcugggggc cacaugagug ug 22
 50
 <210> 577
 <211> 21
 <212> RNA
 <213> Homo sapiens
 55
 <400> 577
 ccaaaucuug aucagaagcc u 21
 <210> 578
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 578
 agggggaugg cagagcaaaa uu 22
 <210> 579
 <211> 21

<212> RNA
 <213> Homo sapiens

5 <400> 579
 auccuugcua ucugggugcu a 21

<210> 580
 <211> 22
 <212> RNA
 10 <213> Homo sapiens

<400> 580
 aaaaguaauu gugguuuuug cc 22

15 <210> 581
 <211> 22
 <212> RNA
 <213> Homo sapiens

20 <400> 581
 aaaaguaauu gcggauuuug cc 22

<210> 582
 <211> 22
 25 <212> RNA
 <213> Homo sapiens

<400> 582
 aaaaguaauu gcggucuuug gu 22

30 <210> 583
 <211> 20
 <212> RNA
 <213> Homo sapiens

35 <400> 583
 cgugccaccc uuuuccccag 20

40 <210> 584
 <211> 17
 <212> RNA
 <213> Homo sapiens

45 <400> 584
 ucccaccgcu gccaccc 17

<210> 585
 <211> 21
 <212> RNA
 50 <213> Homo sapiens

<400> 585
 uggacugccc ugaucuggag a 21

55 <210> 586
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 586 uuuagagacg gggucuugcu cu	22
5	<210> 587 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 587 ucuacagugc acgugucucc ag	22
15	<210> 588 <211> 21 <212> RNA <213> Homo sapiens	
20	<400> 588 agggaggac gggggcugug c	21
25	<210> 589 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 589 ccaguauuaa cugugcugcu ga	22
35	<210> 590 <211> 21 <212> RNA <213> Homo sapiens	
40	<400> 590 ugagugccgg ugccugcccu g	21
45	<210> 591 <211> 21 <212> RNA <213> Homo sapiens	
50	<400> 591 caagucacua gugguuccgu u	21
55	<210> 592 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 592 ucugcaagug ucagaggcga gg	22
	<210> 593 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 593 auuguccuug cuguuaggag au	22

<210> 594
 <211> 22
 <212> RNA
 <213> Homo sapiens
 5
 <400> 594
 ugaccgauuu cuccuggugu uc 22
 10
 <210> 595
 <211> 22
 <212> RNA
 <213> Homo sapiens
 15
 <400> 595
 gauaucagcu caguaggcac cg 22
 20
 <210> 596
 <211> 21
 <212> RNA
 <213> Homo sapiens
 <400> 596
 cacagcaagu guagacaggc a 21
 25
 <210> 597
 <211> 22
 <212> RNA
 <213> Homo sapiens
 30
 <400> 597
 agcuuuuggg aaucaggua gu 22
 35
 <210> 598
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 598
 uauaaaauga gggcaguaag ac 22
 40
 <210> 599
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 599
 aggauuucag aaauacuggu gu 22
 50
 <210> 600
 <211> 22
 <212> RNA
 <213> Homo sapiens
 55
 <400> 600
 cggguggauc acgaugcaau uu 22
 <210> 601
 <211> 22

<212> RNA
 <213> Homo sapiens

5 <400> 601
 ugpgcuggca gggcaagugc ug 22

<210> 602
 <211> 22
 <212> RNA
 10 <213> Homo sapiens

<400> 602
 ugugggaucu ggaggcaucu gg 22

15 <210> 603
 <211> 22
 <212> RNA
 <213> Homo sapiens

20 <400> 603
 uuucuucuua gacauggcaa cg 22

<210> 604
 <211> 22
 25 <212> RNA
 <213> Homo sapiens

<400> 604
 uuagccaaau guccaucuuu ag 22

30 <210> 605
 <211> 23
 <212> RNA
 <213> Homo sapiens

35 <400> 605
 uggagaucca gugcucgccc gau 23

<210> 606
 <211> 22
 40 <212> RNA
 <213> Homo sapiens

<400> 606
 45 ugaaacugga gcgccuggag ga 22

<210> 607
 <211> 21
 <212> RNA
 50 <213> Homo sapiens

<400> 607
 agcggugcuc cugcgggccg a 21

55 <210> 608
 <211> 21
 <212> RNA
 <213> Homo sapiens

	<400> 608 gagguuuggg gaggauuugc u	21
5	<210> 609 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 609 ucaggacacu ucugaacuug ga	22
15	<210> 610 <211> 23 <212> RNA <213> Homo sapiens	
20	<400> 610 uagcagcggg aacaguucug cag	23
25	<210> 611 <211> 24 <212> RNA <213> Homo sapiens	
30	<400> 611 ucugggcaca ggcggaugga cagg	24
35	<210> 612 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 612 caaaaacugc aaauacuuuu gc	22
45	<210> 613 <211> 20 <212> RNA <213> Homo sapiens	
50	<400> 613 aaaacugcag uuacuuuugc	20
55	<210> 614 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 614 ugugucacuc gaugaccacu gu	22
	<210> 615 <211> 21 <212> RNA <213> Homo sapiens	
	<400> 615 acagucugcu gagguuggag c	21

<210> 616
 <211> 23
 <212> RNA
 <213> Homo sapiens
 5
 <400> 616
 ugugcuugcu cgucccgccc gca 23
 10
 <210> 617
 <211> 23
 <212> RNA
 <213> Homo sapiens
 15
 <400> 617
 gggaugguag accggugacg ugc 23
 20
 <210> 618
 <211> 21
 <212> RNA
 <213> Homo sapiens
 <400> 618
 cccaccucc ucucuccu g 21
 25
 <210> 619
 <211> 22
 <212> RNA
 <213> Homo sapiens
 30
 <400> 619
 uggaguguga caaugguguu ug 22
 35
 <210> 620
 <211> 21
 <212> RNA
 <213> Homo sapiens
 <400> 620
 gugggcgggg gcaggugugu g 21
 40
 <210> 621
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 621
 acccgucccg uucgucccg ga 22
 50
 <210> 622
 <211> 22
 <212> RNA
 <213> Homo sapiens
 55
 <400> 622
 cggaugagca aagaaagugg uu 22
 <210> 623
 <211> 22

<212> RNA
 <213> Homo sapiens

5 <400> 623
 cuggacugag ccaugcuacu gg 22

10 <210> 624
 <211> 26
 <212> RNA
 <213> Homo sapiens

15 <400> 624
 gaugaugaug gcagcaaaau cugaaa 26

20 <210> 625
 <211> 22
 <212> RNA
 <213> Homo sapiens

25 <400> 625
 ggcgacaaaa cgagaccug uc 22

30 <210> 626
 <211> 20
 <212> RNA
 <213> Homo sapiens

35 <400> 626
 ucguuugccu uuuucugcuu 20

40 <210> 627
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 627
 ugugagguug gcauuguugu cu 22

50 <210> 628
 <211> 22
 <212> RNA
 <213> Homo sapiens

55 <400> 628
 ccaccucccc ugcaaacguc ca 22

60 <210> 629
 <211> 23
 <212> RNA
 <213> Homo sapiens

65 <400> 629
 uauggcuuuu uauuccuauug uga 23

70 <210> 630
 <211> 23
 <212> RNA
 <213> Homo sapiens

	<400> 630 uuauugcuua agaauacgcg uag	23
5	<210> 631 <211> 21 <212> RNA <213> Homo sapiens	
10	<400> 631 cauaaaguag aaagcacuac u	21
15	<210> 632 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 632 ggugcagugc ugcaucucug gu	22
25	<210> 633 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 633 caggccauau ugugcugccu ca	22
35	<210> 634 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 634 cugccaaauuc cauaggucac ag	22
45	<210> 635 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 635 uaacacuguc ugguaacgau gu	22
55	<210> 636 <211> 21 <212> RNA <213> Homo sapiens	
	<400> 636 gcugggaagg caaagggacg u	21
	<210> 637 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 637 acagcaggca cagacaggca gu	22

<210> 638
 <211> 21
 <212> RNA
 <213> Homo sapiens
 5
 <400> 638
 ucuggcuguu guggugugca a 21
 10
 <210> 639
 <211> 22
 <212> RNA
 <213> Homo sapiens
 15
 <400> 639
 ugccuggaac auaguaggga cu 22
 20
 <210> 640
 <211> 20
 <212> RNA
 <213> Homo sapiens
 <400> 640
 uagaggaagc uguggagaga 20
 25
 <210> 641
 <211> 22
 <212> RNA
 <213> Homo sapiens
 30
 <400> 641
 ucccuucug caggccugcu gg 22
 35
 <210> 642
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 642
 uucagccagg cuagugcagu cu 22
 40
 <210> 643
 <211> 22
 <212> RNA
 <213> Homo sapiens
 45
 <400> 643
 agaaggggug aaauuaaac gu 22
 50
 <210> 644
 <211> 19
 <212> RNA
 <213> Homo sapiens
 55
 <400> 644
 aucgggccc ucgcgccgg 19
 <210> 645
 <211> 23

<212> RNA
 <213> Homo sapiens

5 <400> 645
 ccugggcagc guguggcuga agg 23

10 <210> 646
 <211> 20
 <212> RNA
 <213> Homo sapiens

15 <400> 646
 ucuggccagc uacgucccca 20

20 <210> 647
 <211> 22
 <212> RNA
 <213> Homo sapiens

25 <400> 647
 guggaguccu ggggaugga ga 22

30 <210> 648
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 648
 aaaagcuggg uugagagggc aa 22

40 <210> 649
 <211> 23
 <212> RNA
 <213> Homo sapiens

45 <400> 649
 uuuaguguga uaauggcguu uga 23

50 <210> 650
 <211> 22
 <212> RNA
 <213> Homo sapiens

55 <400> 650
 ucagcaggca ggcuggugca gc 22

60 <210> 651
 <211> 21
 <212> RNA
 <213> Homo sapiens

65 <400> 651
 ugaguguugu cuacgagggc a 21

70 <210> 652
 <211> 22
 <212> RNA
 <213> Homo sapiens

	<400> 652 auuguagaac cuaagauugg cc	22
5	<210> 653 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 653 cuucccccca gaaucuuca uc	22
15	<210> 654 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 654 uuuguucguu cggcucgcgu ga	22
25	<210> 655 <211> 19 <212> RNA <213> Homo sapiens	
30	<400> 655 acuggacuug gaggcagaa	19
35	<210> 656 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 656 gagcaaugua gguagacugu uu	22
45	<210> 657 <211> 22 <212> RNA <213> Homo sapiens	
50	<400> 657 uguguggauc cuggaggagg ca	22
55	<210> 658 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 658 uuugggacug aucuugaugu cu	22
	<210> 659 <211> 20 <212> RNA <213> Homo sapiens	
	<400> 659 gcucggacug agcagguggg	20

	<210> 660	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 660	
	uucgggcugg ccugcugcuc cgg	23
10		
	<210> 661	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 661	
	uaauacuguc ugguaaaacc gu	22
20		
	<210> 662	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 662	
	accugucugu ggaaaggagc ua	22
25		
	<210> 663	
	<211> 17	
	<212> RNA	
	<213> Homo sapiens	
30		
	<400> 663	
	uuggaggcgu ggguuuu	17
35		
	<210> 664	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
40		
	<400> 664	
	ccaggaggcg gaggaggugg ag	22
45		
	<210> 665	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 665	
	cuagugcucu ccguuacaag ua	22
50		
	<210> 666	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 666	
	cgcgcgccg ugcucggagc ag	22
	<210> 667	
	<211> 22	

	<212> RNA	
	<213> Homo sapiens	
5	<400> 667 ugugacaaua gagaugaaca ug	22
	<210> 668	
	<211> 22	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 668 uuuggacaga aaacacgcag gu	22
15	<210> 669	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 669 aacuguuugc agaggaaacu ga	22
	<210> 670	
	<211> 22	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 670 aacucguguu caaagccuuu ag	22
30	<210> 671	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 671 uuucuucuua gacauggcag cu	22
	<210> 672	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 672 cuuccggucu gugagccccg uc	22
	<210> 673	
	<211> 23	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 673 cugggggacg cgugagcgcg age	23
55	<210> 674	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	

	<400> 674 auacauguca gauuguaugc c	21
5	<210> 675 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 675 aacgggaug caggcuguau cu	22
15	<210> 676 <211> 23 <212> RNA <213> Homo sapiens	
20	<400> 676 ucugaguucc uggagccugg ucu	23
25	<210> 677 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 677 gagcaggcga ggcugggcug aa	22
35	<210> 678 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 678 agaagauugc agaguaaguu cc	22
45	<210> 679 <211> 25 <212> RNA <213> Homo sapiens	
50	<400> 679 agcggggagg aagugggcgc ugcuu	25
55	<210> 680 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 680 ugagggcucc aggugacggu gg	22
	<210> 681 <211> 25 <212> RNA <213> Homo sapiens	
	<400> 681 caugcugacc ucccuccugc cccag	25

	<210> 682	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	
5		
	<400> 682	
	ucaggcaaag ggauauuuac aga	23
10		
	<210> 683	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
15		
	<400> 683	
	aaggcccggg cuuuccucc ag	22
20		
	<210> 684	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
25		
	<400> 684	
	ugcggggaca ggccaggga uc	22
30		
	<210> 685	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
35		
	<400> 685	
	agccaggcuc ugaaggga gu	22
40		
	<210> 686	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
45		
	<400> 686	
	cgccugccca gccuccugc u	21
50		
	<210> 687	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
55		
	<400> 687	
	auagcaauug cucuuugga a	21
60		
	<210> 688	
	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
65		
	<400> 688	
	aauguuggaa uccucguag ag	22
70		
	<210> 689	
	<211> 23	

	<212> RNA	
	<213> Homo sapiens	
5	<400> 689 ucugcacugu gaguuggcug gcu	23
	<210> 690	
	<211> 21	
10	<212> RNA	
	<213> Homo sapiens	
	<400> 690 uugaaaggcu auuucuuggu c	21
15	<210> 691	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
20	<400> 691 ugcuguauug ucagguagug a	21
	<210> 692	
	<211> 24	
25	<212> RNA	
	<213> Homo sapiens	
	<400> 692 agggcuggac ucagcggcgg agcu	24
30	<210> 693	
	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
35	<400> 693 aauuugguuu cugaggcacu uagu	24
	<210> 694	
40	<211> 22	
	<212> RNA	
	<213> Homo sapiens	
	<400> 694 ugaggacagg gcaauucac ga	22
	<210> 695	
	<211> 21	
	<212> RNA	
50	<213> Homo sapiens	
	<400> 695 uuucccuuc cauccuggca g	21
55	<210> 696	
	<211> 23	
	<212> RNA	
	<213> Homo sapiens	

	<400> 696 cagggcucag ggauuggaug gag	23
5	<210> 697 <211> 22 <212> RNA <213> Homo sapiens	
10	<400> 697 auucugcauu uuuagcaagu uc	22
15	<210> 698 <211> 22 <212> RNA <213> Homo sapiens	
20	<400> 698 aacggcaaug acuuuuguac ca	22
25	<210> 699 <211> 22 <212> RNA <213> Homo sapiens	
30	<400> 699 gaaaguaauu gcuguuuuug cc	22
35	<210> 700 <211> 22 <212> RNA <213> Homo sapiens	
40	<400> 700 aaaaguauu gcgguuuuug cu	22
45	<210> 701 <211> 21 <212> RNA <213> Homo sapiens	
50	<400> 701 aaaaguauu gcgguuuuug c	21
55	<210> 702 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 702 caagaaccuc aguugcuuuu gu	22
	<210> 703 <211> 22 <212> RNA <213> Homo sapiens	
	<400> 703 auauuaccuu uagcucaucu uu	22

<210> 704
 <211> 22
 <212> RNA
 <213> Homo sapiens
 5
 <400> 704
 uaaaacuuua agugugccua gg 22
 10
 <210> 705
 <211> 22
 <212> RNA
 <213> Homo sapiens
 15
 <400> 705
 cagagugaca agcugguuaa ag 22
 20
 <210> 706
 <211> 21
 <212> RNA
 <213> Homo sapiens
 <400> 706
 acuggcauua gugggacuuu u 21
 25
 <210> 707
 <211> 24
 <212> RNA
 <213> Homo sapiens
 30
 <400> 707
 uacagaugca gauucucuga cuuc 24
 35
 <210> 708
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 708
 cucauuuuag uagucugaug cc 22
 40
 <210> 709
 <211> 19
 <212> RNA
 <213> Homo sapiens
 45
 <400> 709
 uuauugucac guucugauu 19
 50
 <210> 710
 <211> 22
 <212> RNA
 <213> Homo sapiens
 55
 <400> 710
 uucuggauaa caugcugaag cu 22
 <210> 711
 <211> 22

<212> RNA
 <213> Homo sapiens

5 <400> 711
 uugugucaau augcgaugau gu 22

10 <210> 712
 <211> 22
 <212> RNA
 <213> Homo sapiens

15 <400> 712
 cacacacugc aauuacuuuu gc 22

20 <210> 713
 <211> 21
 <212> RNA
 <213> Homo sapiens

25 <400> 713
 cacaagguau ugguauuacc u 21

30 <210> 714
 <211> 22
 <212> RNA
 <213> Homo sapiens

35 <400> 714
 uccauuacac uaccugccu cu 22

40 <210> 715
 <211> 22
 <212> RNA
 <213> Homo sapiens

45 <400> 715
 ugugcgagg gagaccucuc cc 22

50 <210> 716
 <211> 17
 <212> RNA
 <213> Homo sapiens

55 <400> 716
 cgcgccgggc ccggguu 17

50 <210> 717
 <211> 18
 <212> RNA
 <213> Homo sapiens

55 <400> 717
 cggggcggca ggggccuc 18

55 <210> 718
 <211> 18
 <212> RNA
 <213> Homo sapiens

	<400> 718	
	gaggcugaag gaagaugg	18
5	<210> 719	
	<211> 17	
	<212> RNA	
	<213> Homo sapiens	
10	<400> 719	
	gcugggcgag gcuggca	17
	<210> 720	
	<211> 22	
	<212> RNA	
15	<213> Homo sapiens	
	<400> 720	
	guggguacgg cccagugggg gg	22
20	<210> 721	
	<211> 17	
	<212> RNA	
	<213> Homo sapiens	
25	<400> 721	
	gcuucuguag uguaguc	17
	<210> 722	
	<211> 18	
30	<212> RNA	
	<213> Homo sapiens	
	<400> 722	
	cggcuggagg ugugagga	18
35	<210> 723	
	<211> 21	
	<212> RNA	
	<213> Homo sapiens	
40	<400> 723	
	uuacaggcgu gaaccaccgc g	21
	<210> 724	
45	<211> 24	
	<212> RNA	
	<213> Homo sapiens	
	<400> 724	
50	ugggcuaagg gagaugauug ggua	24
	<210> 725	
	<211> 21	
	<212> RNA	
55	<213> Homo sapiens	
	<400> 725	
	ccgucgccgc caccgagcc g	21

<210> 726
 <211> 22
 <212> RNA
 <213> Homo sapiens
 5
 <400> 726
 gguccagagg ggagauaggu uc 22
 10
 <210> 727
 <211> 24
 <212> RNA
 <213> Homo sapiens
 15
 <400> 727
 cauagcccgg ugcugguac auga 24
 20
 <210> 728
 <211> 17
 <212> RNA
 <213> Homo sapiens
 <400> 728
 ccccgccacc gccuugg 17
 25
 <210> 729
 <211> 24
 <212> RNA
 <213> Homo sapiens
 30
 <400> 729
 ugucagugac uccugcccu uggu 24
 35
 <210> 730
 <211> 22
 <212> RNA
 <213> Homo sapiens
 <400> 730
 ucuggggaug aggacagugu gu 22
 40
 <210> 731
 <211> 21
 <212> RNA
 <213> Homo sapiens
 45
 <400> 731
 ucucaguaag uggcacucug u 21
 50
 <210> 732
 <211> 18
 <212> RNA
 <213> Homo sapiens
 55
 <400> 732
 agagcagaag gaugagau 18
 <210> 733
 <211> 22

<212> RNA
 <213> Homo sapiens

 <400> 733
 5 ucagggcagg cacaguggcu ca 22

 <210> 734
 <211> 22
 <212> RNA
 10 <213> Homo sapiens

 <400> 734
 acuuuccuca cucccgugaa gu 22

 15 <210> 735
 <211> 21
 <212> RNA
 <213> Homo sapiens

 20 <400> 735
 ugacagcgcc cugccuggcu c 21

 <210> 736
 <211> 17
 25 <212> RNA
 <213> Homo sapiens

 <400> 736
 aagacugaga ggagggga 17
 30
 <210> 737
 <211> 23
 <212> RNA
 <213> Homo sapiens
 35
 <400> 737
 agaacucuug cagucuuga ugu 23

 <210> 738
 40 <211> 57
 <212> DNA
 <213> Homo sapiens

 <220>
 45 <223> genomic sequence of miRNA precursor

 <400> 738
 ggccgcgccc cgtttccag gacaaagggc actccgcacc ggaccctgtt cccagcg 57
 50
 <210> 739
 <211> 64
 <212> DNA
 <213> Homo sapiens

 55 <220>
 <223> genomic sequence of miRNA precursor

 <400> 739

cccactccct ggcgccgctt gtggagggcc caagtccttc tgattgaggc ccaacccgtg 60
gaag 64

5 <210> 740
<211> 56
<212> DNA
<213> Homo sapiens

10 <220>
<223> genomic sequence of miRNA precursor

<400> 740
cgccgggacc ggggtccggg gcggagtgcc cttctcctg ggaaacgggg tgcggc 56

15 <210> 741
<211> 81
<212> DNA
<213> Homo sapiens

20 <220>
<223> genomic sequence of miRNA precursor

<400> 741

25 gcttcacgtc cccaccggcg gcggcgggcg tggcagtggc ggccggcgcg gcggtggcgg 60
cggcgggcggc ggcggcggt c 81

30 <210> 742
<211> 89
<212> DNA
<213> Homo sapiens

35 <220>
<223> genomic sequence of miRNA precursor

<400> 742

40 gccgcccccg ccgcccgcgc cgcgcgcgcc gccgcgcgcg ccgcccagctt taggatcagg 60
cctcaggtga gtcggagggg ccgggccc 89

45 <210> 743
<211> 22
<212> RNA
<213> Homo sapiens

<400> 743

50 gccggagugc ccuucuuccu gg 22

<210> 744
<211> 22
<212> RNA

55 <213> Homo sapiens

<400> 744
ggagggccca aguccuucug au 22

<210> 745
 <211> 22
 <212> RNA
 <213> Homo sapiens
 5
 <400> 745
 gaccaggguc cggugcggag ug 22
 <210> 746
 10 <211> 17
 <212> DNA
 <213> Homo sapiens
 <220>
 15 <223> genomic primer sequence
 <400> 746
 aggtccggt gcggagt 17
 20 <210> 747
 <211> 12
 <212> RNA
 <213> Bos taurus
 25 <400> 747
 ggguccggug cg 12
 <210> 748
 <211> 20
 30 <212> DNA
 <213> Homo sapiens
 <220>
 <223> genomic primer sequence
 35 <400> 748
 tgcggagtgc ctttgcct 20
 <210> 749
 40 <211> 22
 <212> DNA
 <213> Homo sapiens
 <220>
 45 <223> genomic primer sequence
 <400> 749
 ggagggccca agtccttctg at 22
 50 <210> 750
 <211> 14
 <212> RNA
 <213> Caenorhabditis remanei
 55 <400> 750
 cccaagugcu ucug 14
 <210> 751

<211> 18
 <212> DNA
 <213> Homo sapiens

<220>
 <223> genomic primer sequence

<400> 751
 cggagtgtccc ttcttct 18

<210> 752
 <211> 12
 <212> RNA
 <213> Zea mays

<400> 752
 gugcccuucu uc 12

Claims

1. A neural stem cell exosome for use in a method of treating glioma.
2. The neural stem cell exosome for use according to claim 1, wherein the glioma is glioblastoma.
3. The neural stem cell exosome for use according to claim 2, wherein the exosome (i) inhibits glioblastoma cell migration and/or (ii) induces differentiation of a glioblastoma cell.
4. The neural stem cell exosome for use according to claim 2 or claim 3, wherein the exosome:
 - inhibits glioblastoma cell migration as determined using a transmembrane cell migration assay; and/or
 - induces differentiation as determined by a reduction in nestin expression.
5. The neural stem cell exosome for use according to any one of claims 1-4, wherein the exosome is derived from a neural stem cell that:
 - (a) is proliferating;
 - (b) does not express DCX or GFAP; and/or
 - (c) has been cultured in a multi-compartment bioreactor for less than 4 weeks and optionally no more than 1 week.
6. The neural stem cell exosome for use according to any one of claims 1-5, wherein the exosome is derived from a neural stem cell line, optionally wherein the neural stem cell line is conditionally-immortalised and/or grown in serum free medium, and optionally wherein the neural stem cell line is CTX0E03 having ECACC Accession No. 04091601, STR0C05 having ECACC Accession No.04110301 and HPC0A07 having ECACC Accession No.04092302.
7. The neural stem cell exosome for use according to any preceding claim, wherein the exosome has:
 - (a) a size of between 30 nm and 1000 nm, or between 30 and 200 nm, or between 30 and 100 nm, as determined by electron microscopy; or
 - (b) a density in sucrose of 1.1-1.2 g/ml.
8. The neural stem cell exosome for use according to any preceding claim, comprising RNA, optionally wherein the RNA is mRNA and/or miRNA, optionally wherein the exosome comprises:
 - one, two, three or four of hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 and hsa-miR-4532;
 - one, two, three, four or five of hsa-miR-181a-5p, hsa-miR-1246, hsa-miR-127-3p, hsa-miR-21-5p, and hsa-miR-100-5p;
 - one, two, three, four or five of hsa-miR-181a-5p, hsa-let-7a-5p, hsa-let-7f-5p, hsa-miR-92b-3p, and hsa-miR-

9-5p; or
hsa-miR-486-5p.

9. The neural stem cell exosome for use according to any preceding claim, comprising at least 10 of the proteins present in Table 20 or Table 22; and/or comprising one or more of:

(a) a lipid selected from ceramide, cholesterol, sphingomyelin, phosphatidylserine, phosphatidylinositol, and/or phosphatidylcholine;

(b) miRNA, optionally selected from hsa-let-7g, hsa-miR-101, hsa-miR-10a, hsa-miR-10b, hsa-miR-126, hsa-miR-128, hsa-miR-129-5p, hsa-miR-130a, hsa-miR-134, hsa-miR-137, hsa-miR-155, hsa-miR-15a, hsa-miR-15b, hsa-miR-16, hsa-miR-17, hsa-miR-182, hsa-miR-183, hsa-miR-185, hsa-miR-18b, hsa-miR-192, hsa-miR-194, hsa-miR-195, hsa-miR-20a, hsa-miR-20b, hsa-miR-210, hsa-miR-218, hsa-miR-301a, hsa-miR-302a, hsa-miR-302c, hsa-miR-345, hsa-miR-375, hsa-miR-378, hsa-miR-7, hsa-miR-9, hsa-miR-93, hsa-miR-96, and hsa-miR-99a;

(c) a tetraspanin, optionally selected from CD63, CD81, CD9, CD53, CD82 and/or CD37;

(d) TSG101, Alix, CD109 and/or thy-1; and/or

(e) CD133.

10. The neural stem cell exosome for use according to any of claims 1 to 9, wherein the glioma is treated by inhibiting migration of the cancer cells or wherein the glioma is treated by inducing differentiation of the cancer cells; optionally wherein the exosome is as defined in claim 5

Patentansprüche

1. Neuronales Stammzellenexosom zur Anwendung in einem Verfahren zur Behandlung von Gliomen.

2. Neuronales Stammzellenexosom zur Anwendung nach Anspruch 1, wobei das Gliom ein Glioblastom ist.

3. Neuronales Stammzellenexosom zur Anwendung nach Anspruch 2, wobei das Exosom (i) Zellmigration von Glioblastomen hemmt und/oder (ii) die Differenzierung einer Glioblastomzelle bewirkt.

4. Neuronales Stammzellenexosom zur Anwendung nach einem der Ansprüche 2 oder 3, wobei das Exosom:

Zellmigration von Glioblastomen wie vorgegeben unter Anwendung eines transmembranen Zellmigrationstests hemmt; und/oder
die Differenzierung wie vorgegeben durch eine Nestinexpressionreduzierung bewirkt.

5. Neuronales Stammzellenexosom zur Anwendung nach einem der Ansprüche 1-4, wobei das Exosom von einer neuronalen Stammzelle abgeleitet ist, die:

(a) sich vermehrt;

(b) weder DCX noch GFAP exprimiert; und/oder

(c) weniger als 4 Wochen und gegebenenfalls nicht mehr als 1 Woche in einem Bioreaktor mit mehreren Abteilungen kultiviert wurde.

6. Neuronales Stammzellenexosom zur Anwendung nach einem der Ansprüche 1-5, wobei das Exosom von einer neuronalen Stammzellenlinie abgeleitet ist, wobei die neuronale Stammzellenlinie gegebenenfalls bedingt verewigt und/oder in einem serumfreien Medium gewachsen ist, und wobei die neuronale Stammzellenlinie gegebenenfalls CTX0E03 mit der ECACC-Zugangsnummer 04091601, STR0C05 mit der ECACC-Zugangsnummer 04110301 und HPC0A07 mit der ECACC-Zugangsnummer 04092302 ist.

7. Neuronales Stammzellenexosom zur Anwendung nach einem der vorhergehenden Ansprüche, wobei das Exosom:

(a) eine Größe zwischen 30 nm und 1000 nm oder zwischen 30 und 200 nm oder zwischen 30 und 100 nm, wie durch Elektronenmikroskopie vorgegeben, aufweist; oder

(b) eine Saccharose-Dichte von 1,1-1,2 g/ml aufweist.

8. Neuronales Stammzellenexosom zur Anwendung nach einem der vorhergehenden Ansprüche, das RNA umfasst, wobei das RNA gegebenenfalls mRNA und/oder miRNA ist, wobei das Exosom gegebenenfalls Folgendes umfasst:

ein, zwei, drei oder vier hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 und hsa-miR-4532;
 ein, zwei, drei, vier oder fünf hsa-miR-181a-5p, hsa-miR-1246, hsa-miR-127-3p, hsa-miR-21-5p, und hsa-miR-100-5p;
 ein, zwei, drei, vier oder fünf hsa-miR-181a-5p, hsa-let-7a-5p, hsa-let-7f-5p, hsa-miR-92b-3p, und hsa-miR-9-5p; oder hsa-miR-486-5p.

9. Neuronales Stammzellenexosom zur Anwendung nach einem der vorhergehenden Ansprüche, das mindestens 10 der in Tabelle 20 oder Tabelle 22 angegebenen Proteine umfasst; und/oder eines oder mehrere der Folgenden umfasst:

(a) ein aus Ceramid, Cholesterin, Sphingomyelin, Phosphatidylserin, Phosphatidylinositol und/oder Phosphatidylcholin ausgewähltes Lipid;
 (b) miRNA, gegebenenfalls ausgewählt aus hsa-let-7g, hsa-miR-101, hsa-miR-10a, hsa-miR-10b, hsa-miR-126, hsa-miR-128, hsa-miR-129-5p, hsa-miR-130a, hsa-miR-134, hsa-miR-137, hsa-miR-155, hsa-miR-15a, hsa-miR-15b, hsa-miR-16, hsa-miR-17, hsa-miR-182, hsa-miR-183, hsa-miR-185, hsa-miR-18b, hsa-miR-192, hsa-miR-194, hsa-miR-195, hsa-miR-20a, hsa-miR-20b, hsa-miR-210, hsa-miR-218, hsa-miR-301a, hsa-miR-302a, hsa-miR-302c, hsa-miR-345, hsa-miR-375, hsa-miR-378, hsa-miR-7, hsa-miR-9, hsa-miR-93, hsa-miR-96, und hsa-miR-99a;
 (c) ein Tetraspanin, das gegebenenfalls aus CD63, CD81, CD9, CD53, CD82 und/oder CD37 ausgewählt ist;
 (d) TSG101, Alix, CD109 und/oder thy-1; und/oder
 (e) CD133.

10. Neuronales Stammzellenexosom zur Anwendung nach einem der Ansprüche 1-9, wobei das Gliom durch Hemmen der Migration von Krebszellen behandelt wird oder wobei das Gliom durch Bewirken einer Differenzierung der Krebszellen behandelt wird; wobei das Exosom gegebenenfalls nach Anspruch 5 definiert wird.

Revendications

1. Exosome de cellules souches neurales destiné à être utilisé dans un procédé de traitement du gliome.
2. Exosome de cellules souches neurales destiné à être utilisé selon la revendication 1, dans lequel le gliome est le glioblastome.
3. Exosome de cellules souches neurales destiné à être utilisé selon la revendication 2, dans lequel l'exosome (i) inhibe la migration de cellules de glioblastome et/ou (ii) induit la différenciation d'une cellule de glioblastome.
4. Exosome de cellules souches neurales destiné à être utilisé selon la revendication 2 ou la revendication 3, dans lequel l'exosome :
 inhibe la migration de cellules de glioblastome telle que déterminée en utilisant un essai de migration de cellules transmembranaires ; et/ou
 induit la différenciation telle que déterminée par une réduction en termes d'expression de nestine.
5. Exosome de cellules souches neurales destiné à être utilisé selon l'une quelconque des revendications 1 à 4, dans lequel l'exosome est dérivé d'une cellule souche neurale qui :
 (a) prolifère ;
 (b) n'exprime pas le DCX ou le GFAP ; et/ou
 (c) a été mise en culture dans un bioréacteur à compartiments multiples pendant moins de 4 semaines et éventuellement pas plus de 1 semaine.
6. Exosome de cellules souches neurales destiné à être utilisé selon l'une quelconque des revendications 1 à 5, dans lequel l'exosome est dérivé d'une lignée de cellules souches neurales, dans lequel la lignée de cellules souches neurales est éventuellement immortalisée sous condition et/ou cultivée dans un milieu exempt de sérum, et dans

lequel la lignée de cellules souches neurales est éventuellement CTX0E03 ayant le numéro d'accès 04091601 à l'ECACC, STR0C05 ayant le numéro d'accès 04110301 à l'ECACC et HPC0A07 ayant le numéro d'accès 04092302 à l'ECACC.

- 5 7. Exosome de cellules souches neurales destiné à être utilisé selon l'une quelconque des revendications précédentes, dans lequel l'exosome comporte :

(a) une taille comprise entre 30 nm et 1 000 nm, ou entre 30 et 200 nm, ou entre 30 et 100 nm, telle que déterminée par microscopie électronique ; ou

10 (b) une densité dans le sucrose comprise entre 1,1 et 1,2 g/ml.

8. Exosome de cellules souches neurales destiné à être utilisé selon l'une quelconque des revendications précédentes, comprenant de l'ARN, dans lequel l'ARN est éventuellement l'ARNm et/ou l'ARNmi, dans lequel l'exosome comprend éventuellement :

15 un, deux, trois ou quatre parmi hsa-miR-1246, hsa-miR-4492, hsa-miR-4488 et hsa-miR-4532 ;

un, deux, trois, quatre ou cinq parmi hsa-miR-181a-5p, hsa-miR-1246, hsa-miR-127-3p, hsa-miR-21-5p, et hsa-miR-100-5p ;

20 un, deux, trois, quatre ou cinq parmi hsa-miR-181a-5p, hsa-let-7a-5p, hsa-let-7f-5p, hsa-miR-92b-3p, et hsa-miR-9-5p ; ou hsa-miR-486-5p.

9. Exosome de cellules souches neurales destiné à être utilisé selon l'une quelconque des revendications précédentes, comprenant au moins 10 des protéines présentes dans le Tableau 20 ou dans le Tableau 22 ; et/ou comprenant l'au moins un parmi :

25 (a) un lipide choisi parmi le céramide, le cholestérol, la sphingomyéline, la phosphatidylsérine, le phosphatidylinositol, et/ou la phosphatidylcholine ;

(b) l'ARNmi, éventuellement choisi parmi hsa-let-7g, hsa-miR-101, hsa-miR-10a, hsa-miR-10b, hsa-miR-126, hsa-miR-128, hsa-miR-129-5p, hsa-miR-130a, hsa-miR-134, hsa-miR-137, hsa-miR-155, hsa-miR-15a, hsa-miR-15b, hsa-miR-16, hsa-miR-17, hsa-miR-182, hsa-miR-183, hsa-miR-185, hsa-miR-18b, hsa-miR-192, hsa-miR-194, hsa-miR-195, hsa-miR-20a, hsa-miR-20b, hsa-miR-210, hsa-miR-218, hsa-miR-301a, hsa-miR-302a, hsa-miR-302c, hsa-miR-345, hsa-miR-375, hsa-miR-378, hsa-miR-7, hsa-miR-9, hsa-miR-93, hsa-miR-96, et hsa-miR-99a ;

(c) une tétraspanine, éventuellement choisie parmi CD63, CD81, CD9, CD53, CD82 et/ou CD37 ;

35 (d) TSG101, Alix, CD109 et/ou thy-1 ; et/ou

(e) CD133.

10. Exosome de cellules souches neurales destiné à être utilisé selon l'une quelconque des revendications 1 à 9, dans lequel le gliome est traité en inhibant la migration des cellules cancéreuses ou dans lequel le gliome est traité en induisant la différenciation des cellules cancéreuses ; dans lequel l'exosome est éventuellement tel que défini selon la revendication 5.







Figure 1







A= CTX0E03 conditioned media
 B= 2ug/ml purified exosomes
 C= Control

Figure 3A






Figure 3B

Figure 3C






Figure 3D




Figure 4




Figure 5

Figure 6




Figure 7A





Figure 7B






Figure 8A




Figure 8B




Figure 9




Figure 10A




Figure 10B




Figure 11

Fold-Regulation of miRNA in Microparticles

▨ CTX0E03 STD ■ CTX0E03 3 week integra






Figure 12

CTX0E03 07EH READ.COUNTS



MVCTX0E03 07EH READ.COUNTS



EXO CTX0E03 07EH READ.COUNTS






Figure 13A

CTX0E03 07EI READ.COUNTS



MVCTX0E03 07EI READ.COUNTS



EXO CTX0E03 07EI READ.COUNTS






Figure 13B

Figure 13C


EXO 07EH




EXO 6W



EXO1 W11



EXO2 W11



EXO3 W11




Figure 13D – Standard Culture




Figure 13E – Exo 6 Week




Figure 13F




Figure 13G





Figure 13H






Figure 14




Figure 15




Figure 16

Figure 17A




Figure 17B




Figure 17C




Figure 17D




Figure 18




Figure 19




Figure 20





Figure 21






Figure 22A




Figure 22B







Figure 22C



U373MG 20nM miRNA 2500 cells/well 10%FBS




Figure 23B




Figure 23C

U87 20nM miRNA 9000cells/well basal




Figure 24: Tumour volume at assignment




Figure 25: Mean body weights





Figure 26: Mean tumour volume by study day¹

Figure 27: Mean tumour volume by study day – Day 25 truncated format





Figure 28: Final tumour weights






Figure 29: Survival plot for the dosing phase

Figure 30: Individual body weights – Absolute

Absolute				Date/Study Day																					
				Pre-Dosing		Dosing & Recovery																			
						4/30	5/2	5/6	5/7	5/8	5/9	5/10	5/11	5/12	5/13	5/14	5/15	5/16	5/17	5/18	5/19	5/20	5/21	5/22	5/23
Animal ID	Group	Sex	Fate	Code	Day	6	8	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29
00077E6CC9-7	1	F	TS	24.3	24.3	25.4	25.7	26.4	26.4	26.4	26.4	25.5	25.5	26.4	26.4	27.0	27.0	27.0	28.0	28.3					
00077E7D9C-3	1	F	TS	25	22.1	23.6	23.3	23.3	23.3	23.3	23.3	23.5	24.1	24.4	25.4	25.4	26.2	26.1	26.1	26.1	26.6	26.6	26.6	26.6	
00077E7ED6-17	1	F	TS	41	22.7	24.0	25.0	25.5	25.2	24.0	25.2	24.1	25.2	26.0	26.7	26.8	27.7	26.9	26.7	27.3	27.2	27.2	27.1	27.1	
00077E7FDB-14	1	F	SS	50	26.7	26.1	26.7	26.8	26.8	26.8	26.8	26.0	27.7	27.7	26.9	26.9	26.7	27.3	27.3	27.2	27.2	27.2	27.1	27.1	
00077E8B08-20	1	F	TS	27	22.3	22.5	23.7	23.9	23.4	23.4	23.4	24.1	24.8	24.8	24.8	24.8	25.1	24.8	25.1	24.7	24.7	24.7	24.7	24.7	
00077E6AFE-21	2	F	SS	50	20.8	20.8	21.4	22.1	22.3	22.3	22.3	21.6	22.1	22.1	22.5	22.5	22.5	22.7	22.4	22.4	21.9	21.9	21.9	21.9	
00077E6C4D-12	2	F	SS	50	22.0	22.6	23.3	23.8	24.3	24.3	24.3	24.3	24.7	24.9	24.9	25.3	24.9	25.1	25.1	25.1	25.1	25.1	25.1	25.1	
00077E7EF9-18	2	F	TS	32	24.0	24.1	24.9	25.1	24.8	23.9	24.8	23.9	25.0	25.0	24.8	26.2	26.5	26.2	26.5	27.5	27.5	27.1	27.1	27.1	
00077E8B5F-30	2	F	TS	25	23.5	23.2	25.3	25.5	25.2	25.2	25.2	25.2	26.1	27.3	27.3	27.3	27.5	27.5	27.3						
00077E8C0D-4	2	F	TS	22	22.9	23.3	24.6	24.4	24.6	23.9	23.9	23.9	24.6	24.6	25.3	25.3	25.3	25.3	25.3						
00077E6ABD-11	3	F	TS	25	23.0	22.8	23.9	23.8	24.0	23.2	23.2	23.2	24.6	24.6	25.5	25.5	25.5	25.5	25.1	25.1	27.1	27.1	27.1	27.1	
00077E6B0F-22	3	F	TS	27	24.4	24.0	25.1	25.7	25.9	25.9	25.9	24.8	25.9	26.9	27.3	27.3	27.2	27.2	28.3	28.3	28.3	28.3	28.3	28.3	
00077E6D43-16	3	F	TS	27	26.7	26.1	26.6	26.5	26.5	26.5	26.5	25.4	26.9	27.3	27.3	27.1	27.5	27.5	27.6	28.0	28.0	28.0	28.0	28.0	
00077E7DF2-13	3	F	TS	27	24.1	24.1	25.3	25.2	25.6	25.6	25.6	25.3	26.9	26.9	27.1	27.1	27.5	27.5	27.6	28.0	28.0	28.0	28.0	28.0	
00077E8C4C-5	3	F	TS	32	21.7	21.4	22.6	22.3	22.9	22.3	22.3	22.3	23.7	23.7	23.7	24.2	24.2	24.4	24.4	24.6	24.6	24.6	24.6	24.6	
00077E6B73-26	4	F	TS	25	24.4	23.5	24.8	24.2	25.0	24.1	24.1	24.1	25.8	25.8	26.0	26.6	26.7	26.6	26.7	26.7	26.7	26.7	26.7	26.7	
00077E6D2A-28	4	F	TS	27	22.5	22.1	23.3	23.0	23.2	23.2	23.1	23.1	23.9	23.9	23.5	23.9	24.1	24.1	25.1	25.1	25.1	25.1	25.1	25.1	
00077E6D5A-2	4	F	TS	34	23.2	25.7	24.6	24.0	23.8	24.6	24.6	23.7	25.0	25.0	24.2	24.7	24.6	24.7	24.6	24.8	24.8	24.8	24.8	24.8	
00077E6D6E-23	4	F	TS	25	25.5	25.9	26.0	26.6	26.6	26.6	26.6	25.9	27.3	27.3	28.1	28.4	28.5	28.4	28.5	28.5	28.5	28.5	28.5	28.5	
00077E7DCC-29	4	F	TS	39	24.3	24.3	25.8	25.5	25.5	25.5	25.5	24.4	26.5	26.5	26.5	26.3	26.4	26.3	26.4	27.3	27.3	27.1	27.1	27.1	
00077E6B14-1	5	F	SS	50	21.8	21.3	22.5	22.2	22.4	22.2	22.3	21.6	22.3	22.6	22.5	22.7	22.9	23.2	22.5	23.1	23.1	22.9	23.4	23.4	
00077E6B00-27	5	F	TS	41	23.7	24.1	24.6	24.9	24.3	24.6	24.9	24.1	25.0	25.4	25.2	25.8	25.0	25.1	24.8	24.9	24.9	24.7	24.4	24.4	
00077E80AD-19	5	F	TS	48	23.0	22.7	23.4	22.2	23.1	23.0	23.8	22.6	23.4	23.7	24.1	25.1	25.1	25.6	24.6	24.5	24.6	24.3	24.7	24.7	
00077E8A89-25	5	F	TS	22	22.3	22.4	23.1	22.8	23.3	23.5	23.3	23.6	23.0	23.9	24.6	24.8	25.0	24.8	25.0	24.6	24.5	24.6	24.3	24.7	
00077E8C04-9	5	F	SS	50	25.3	25.3	25.6	25.5	26.2	26.0	26.5	25.8	26.8	27.1	27.3	27.3	26.8	27.2	26.5	26.4	26.4	26.2	26.5	26.5	
00077E6C09-7	1	F	TS	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	50	50	
00077E7D9C-3	1	F	TS	27	27.1	28.1	27.1	26.5	26.3	27.8	26.8	26.8	27.8	27.8	27.0	28.2	27.9	28.8	28.8	29.1	29.3	29.3	29.3	29.3	
00077E7ED6-17	1	F	TS	41	27.1	28.1	27.1	26.5	26.3	27.8	26.8	26.8	27.8	27.8	27.0	28.2	27.9	28.8	28.8	29.1	29.3	29.3	29.3	29.3	
00077E7FDB-14	1	F	SS	50	26.7	26.1	26.7	26.8	26.8	26.8	26.8	26.0	27.7	27.7	26.9	26.9	26.7	27.3	27.3	27.2	27.2	27.2	27.1	27.1	
00077E8B08-20	1	F	TS	27	22.3	22.5	23.7	23.9	23.4	23.4	23.4	24.1	24.8	24.8	24.8	24.8	25.1	24.8	25.1	24.7	24.7	24.7	24.7	24.7	
00077E6AFE-21	2	F	SS	50	20.8	20.8	21.4	22.1	22.3	22.3	22.3	21.6	22.1	22.1	22.5	22.5	22.5	22.7	22.4	22.4	21.9	21.9	21.9	21.9	
00077E6C4D-12	2	F	SS	50	22.0	22.6	23.3	23.8	24.3	24.3	24.3	24.3	24.7	24.9	24.9	25.3	24.9	25.1	25.1	25.1	25.1	25.1	25.1	25.1	
00077E7EF9-18	2	F	TS	32	24.0	24.1	24.9	25.1	24.8	23.9	24.8	23.9	25.0	25.0	24.8	26.2	26.5	26.2	26.5	27.5	27.5	27.1	27.1	27.1	
00077E8B5F-30	2	F	TS	25	23.5	23.2	25.3	25.5	25.2	25.2	25.2	25.2	26.1	27.3	27.3	27.3	27.5	27.5	27.3						
00077E8C0D-4	2	F	TS	22	22.9	23.3	24.6	24.4	24.6	23.9	23.9	23.9	24.6	24.6	25.3	25.3	25.3	25.3	25.3						
00077E6ABD-11	3	F	TS	25	23.0	22.8	23.9	23.8	24.0	23.2	23.2	23.2	24.6	24.6	25.5	25.5	25.5	25.5	25.1	25.1	27.1	27.1	27.1	27.1	
00077E6B0F-22	3	F	TS	27	24.4	24.0	25.1	25.7	25.9	25.9	25.9	24.8	25.9	26.9	27.3	27.3	27.2	27.2	28.3	28.3	28.3	28.3	28.3	28.3	
00077E6D43-16	3	F	TS	27	26.7	26.1	26.6	26.5	26.5	26.5	26.5	25.4	26.9	27.3	27.3	27.1	27.5	27.5	27.6	28.0	28.0	28.0	28.0	28.0	
00077E7DF2-13	3	F	TS	27	24.1	24.1	25.3	25.2	25.6	25.6	25.6	25.3	26.9	26.9	27.1	27.1	27.5	27.5	27.6	28.0	28.0	28.0	28.0	28.0	
00077E8C4C-5	3	F	TS	32	21.7	21.4	22.6	22.3	22.9	22.3	22.3	22.3	23.7	23.7	23.7	24.2	24.2	24.4	24.4	24.6	24.6	24.6	24.6	24.6	
00077E6B73-26	4	F	TS	25	24.4	23.5	24.8	24.2	25.0	24.1	24.1	24.1	25.8	25.8	26.0	26.6	26.7	26.6	26.7	26.7	26.7	26.7	26.7	26.7	
00077E6D2A-28	4	F	TS	27	22.5	22.1	23.3	23.0	23.2	23.2	23.1	23.1	23.9	23.9	23.5	23.9	24.1	24.1	25.1	25.1	25.1	25.1	25.1	25.1	
00077E6D5A-2	4	F	TS	34	23.2	25.7	24.6	24.0	23.8	24.6	24.6	23.7	25.0	25.0	24.2	24.7	24.6	24.7	24.6	24.8	24.8	24.8	24.8	24.8	
00077E6D6E-23	4	F	TS	25	25.5	25.9	26.0	26.6	26.6	26.6	26.6	25.9	27.3	27.3	28.1	28.4	28.5	28.4	28.5	28.5	28.5	28.5	28.5	28.5	
00077E7DCC-29	4	F	TS	39	24.3	24.3	25.8	25.5	25.5	25.5	25.5	24.4	26.5	26.5	26.5	26.3	26.4	26.3	26.4	27.3	27.3	27.1	27.1	27.1	
00077E6ABD-11	3	F	TS	25	21.8	21.3	22.5	22.2	22.4	22.2	22.3	21.6	22.3	22.6	22.5	22.7	22.9	23.2	22.5	23.1	23.1	22.9	23.4	23.4	
00077E6B0F-22	3	F	TS	27	24.4	24.0	25.1	25.7	25.9	25.9	25.9	24.8	25.9	26.9	27.3	27.3	27.2	27.2</							

Figure 31: Individual body weights – Relative

Relative			Pre-Dosing			Date/Study Day																			
						Dosing & Recovery																			
Animal ID	Group	Sex	Fate	Day	4/30	5/2	8	5/6	5/7	5/8	5/9	5/10	5/11	5/12	5/13	5/14	5/15	5/16	5/17	5/18	5/19	5/20	5/21	5/22	5/23
00077E6CC9-7	1	F	TS	27	95.7%	95.7%		100.0%	101.2%	103.9%			100.4%			103.9%	106.3%	106.3%	106.3%	106.3%	110.2%		111.4%		
00077E7D9C-3	1	F	TS	25	93.6%	93.6%	100.0%	98.7%	98.7%	98.7%			99.6%			103.4%	107.6%	111.0%	110.6%	107.6%					
00077E7ED6-17	1	F	TS	41	90.8%	96.0%	100.0%	102.0%	100.8%	100.4%			96.4%			99.2%	100.8%	102.8%	103.6%	100.8%	103.6%		106.4%		106.4%
00077E7FDB-14	1	F	SS	50	100.0%	97.8%	100.0%	100.0%	100.4%	100.0%			97.4%			103.7%	100.7%	100.0%	102.2%	100.7%	102.2%		101.9%		101.5%
00077E8B08-20	1	F	TS	27	94.1%	94.9%	100.0%	100.8%	98.7%				101.7%			104.6%	104.6%	105.9%	104.6%	105.9%	104.6%		104.2%		
00077E8AFE-21	2	F	SS	50	97.2%	97.2%	100.0%	103.3%	104.2%				100.9%			103.3%	105.1%	106.1%	106.9%	105.1%	106.1%		104.7%		102.3%
00077E8C4D-12	2	F	SS	50	94.4%	97.0%	100.0%	102.1%	104.3%				104.3%			106.0%	106.9%	106.9%	106.9%	106.9%	107.7%		107.7%		107.7%
00077E7EF9-18	2	F	TS	32	96.4%	96.8%	100.0%	100.8%	99.6%				99.0%			100.4%	99.6%	105.2%	106.4%	105.2%	106.4%		110.4%		108.8%
00077E8B5F-30	2	F	TS	25	92.9%	91.7%	100.0%	100.8%	99.6%				99.6%			103.2%	107.9%	107.9%	107.9%	107.9%	107.9%		108.7%		108.7%
00077E8CUD-4	2	F	TS	22	93.1%	94.7%	100.0%	99.2%	100.0%				97.2%			100.0%	102.8%	102.8%	102.8%	102.8%					
00077E8ABD-11	3	F	TS	25	96.2%	95.4%	100.0%	99.6%	100.4%				97.1%			102.9%	102.9%	106.7%	105.0%	102.9%	105.0%				
00077E8B0F-22	3	F	TS	27	97.2%	95.6%	100.0%	102.4%	103.2%				98.8%			103.2%	104.4%	106.4%	104.8%	104.4%	104.8%		108.0%		
00077E8D43-16	3	F	TS	27	100.4%	98.1%	100.0%	97.4%	99.6%				95.5%			101.1%	102.6%	102.3%	102.3%	102.6%	102.3%		106.4%		106.4%
00077E7DF2-13	3	F	TS	27	95.3%	95.3%	100.0%	99.6%	101.2%				100.0%			106.3%	107.1%	108.7%	109.1%	107.1%	109.1%		110.7%		110.7%
00077E8C4C-5	3	F	TS	32	96.0%	94.7%	100.0%	98.7%	101.3%				98.7%			104.9%	104.9%	107.3%	108.0%	104.9%	108.0%		108.8%		107.1%
00077E8E73-26	4	F	TS	25	98.4%	94.8%	100.0%	97.6%	100.8%				97.2%			104.0%	104.8%	107.3%	107.7%	104.8%	107.7%				
00077E8D2A-28	4	F	TS	27	96.6%	94.8%	100.0%	98.7%	98.6%				99.1%			102.6%	100.9%	102.6%	103.4%	100.9%	103.4%		107.7%		102.8%
00077E6D5A-2	4	F	TS	34	94.3%	104.5%	100.0%	97.6%	96.7%				96.3%			101.6%	98.4%	100.4%	100.0%	98.4%	100.0%		100.8%		
00077E6DEE-23	4	F	TS	25	98.1%	99.6%	100.0%	102.3%	102.3%				99.6%			105.0%	108.1%	109.2%	109.6%	108.1%	109.6%				
00077E7EDC-29	4	F	TS	39	94.2%	94.2%	100.0%	98.8%	98.8%				94.6%			102.7%	102.7%	101.9%	102.3%	102.7%	102.3%		105.8%		105.0%
00077E6B81-14	5	F	SS	50	96.9%	94.7%	100.0%	98.7%	98.2%	99.6%	98.7%	99.1%	96.0%	99.1%	100.4%	100.0%	100.9%	101.8%	103.1%	100.0%	102.7%	102.7%	101.8%	104.0%	
00077E6BD0-27	5	F	TS	41	96.3%	98.0%	100.0%	101.2%	98.8%	100.0%	101.2%	101.2%	98.0%	101.6%	103.3%	102.4%	104.9%	101.6%	102.0%	100.8%	101.2%	101.2%	100.4%	99.2%	
00077E8AD0-19	5	F	TS	48	98.3%	97.0%	100.0%	94.9%	98.7%	98.3%	101.7%	100.0%	96.6%	100.0%	103.0%	107.3%	103.0%	107.3%	109.4%	105.1%	104.7%	105.1%	103.8%	105.6%	
00077E8A89-25	5	F	TS	22	96.5%	97.0%	100.0%	98.7%	100.9%	101.7%	100.9%	102.2%	99.6%	103.5%	106.5%	107.4%	108.2%								
00077E8CC4-9	5	F	SS	50	98.8%	98.8%	100.0%	99.6%	102.3%	101.6%	103.5%	103.9%	100.8%	104.7%	105.9%	106.6%	106.6%	104.7%	106.3%	103.5%	103.1%	103.1%	102.3%	103.5%	
Animal ID	Group	Sex	Fate	Day	5/24	5/25	5/26	5/27	5/28	5/29	5/30	5/31	6/1	6/2	6/3	6/4	6/5	6/6	6/7	6/8	6/9	6/10	6/11	6/13	
00077E6CC9-7	1	F	TS	27	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	50	
00077E7D9C-3	1	F	TS	25			108.4%		106.0%		105.2%			107.2%		108.0%									
00077E7ED6-17	1	F	TS	41								104.1%		104.1%		105.6%									
00077E7FDB-14	1	F	SS	50			105.2%	101.5%													107.9%	109.0%	109.7%		
00077E8B08-20	1	F	TS	27															104.5%						
00077E8AFE-21	2	F	SS	50			103.7%	96.7%			100.5%		101.4%	101.4%		103.3%	104.2%				105.1%	108.9%	106.1%		
00077E8C4D-12	2	F	SS	50			107.3%	103.4%			105.6%		105.6%	105.6%		106.9%	109.0%				106.9%	110.3%	105.2%		
00077E7EF9-18	2	F	TS	32			112.0%																		
00077E8B5F-30	2	F	TS	25																					
00077E8CUD-4	2	F	TS	22			111.1%																		
00077E8ABD-11	3	F	TS	25																					
00077E8B0F-22	3	F	TS	27																					
00077E8D43-16	3	F	TS	27																					
00077E7DF2-13	3	F	TS	27																					
00077E8C4C-5	3	F	TS	32																					
00077E6B873-26	4	F	TS	25																					
00077E6D2A-28	4	F	TS	27																					
00077E6D5A-2	4	F	TS	34																					
00077E6DEE-23	4	F	TS	25																					
00077E7DCC-29	4	F	TS	39																					
00077E6B81-14	5	F	SS	50	105.8%	105.8%	105.3%	104.4%	104.4%	102.2%	107.6%	107.6%	106.2%	102.7%	103.1%	103.6%	105.8%	106.2%	107.6%	109.8%	108.0%	104.0%	105.8%	104.9%	
00077E6BD0-27	5	F	TS	41	98.4%	100.8%	98.4%	97.2%	95.9%	93.5%	96.3%	95.9%	96.3%	93.1%	94.7%	95.1%									
00077E80AD-19	5	F	TS	48	104.7%	106.8%	108.1%	105.6%	107.7%	104.7%	109.0%	111.5%	106.4%	109.8%	108.5%	109.4%	110.7%	110.7%	110.3%	111.5%	108.1%	109.8%			
00077E8A89-25	5	F	TS	22																					
00077E8CC4-9	5	F	SS	50	105.9%	104.7%	105.9%	102.0%	105.1%	102.0%	105.5%	107.4%	107.4%	104.3%	105.5%	107.6%	106.6%	107.4%	108.2%	110.2%	108.6%	103.9%	108.6%		

Figure 32: Individual tumour volume measurements

			Date/Study Day																												
			Pre-Dosing		Dosing & Recovery																										
Animal ID	Group	Sex	Pace	Code	Day	4/30	5/2	5/6	5/7	5/9	5/11	5/13	5/14	5/16	5/19	5/21	5/23	5/26	5/28	5/30	5/32	5/34	5/36	5/38	5/40	5/42	5/44	5/46	5/48	5/50	
000077E0C0-7	1	F	TS	27	65.0	121.1	146.1	140.3	227.3	343.2	571.4	612.1	870.7	1,248.1	1,603.0																
000077E09C-3	1	F	TS	25	152.6	176.5	198.4	233.9	374.7	529.9	674.7	998.4	1,313.2	1,765.0																	
000077E0ED-6	1	F	TS	41	55.2	92.3	106.0	132.6	195.0	221.2	343.2	400.0	537.8	740.7	819.3	1,107.1	1,225.3	1,306.8	1,278.3	1,412.4	1,661.0										
000077E0FD-14	1	F	SS	50	78.7	94.1	181.8	226.5	228.1	216.3	194.5	173.8	122.4	37.0	20.7	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
000077E080-20	1	F	TS	27	55.5	108.8	142.2	172.7	192.8	287.1	372.1	467.0	676.0	1,030.2	1,561.6																
000077E04E-21	2	F	SS	50	66.7	100.0	101.1	131.4	163.3	134.0	124.2	113.7	71.3	36.1	23.7	22.3	17.9	10.6	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
000077E0C0-12	2	F	SS	50	116.6	134.0	248.4	305.9	461.3	509.9	548.4	595.8	619.0	681.2	771.5	777.7	643.5	578.2	497.7	433.0	354.5	363.3	271.5	317.5	261.6						
000077E0FF-18	2	F	TS	32	65.7	121.8	135.0	154.8	198.2	244.5	391.6	451.4	624.3	938.5	1,175.8	1,376.3	1,679.8														
000077E085-30	2	F	TS	25	111.4	157.7	183.0	216.6	295.7	590.6	831.7	916.2	1,168.9	1,687.0																	
000077E0C0-4	2	F	TS	22	321.8	131.2	153.2	219.6	284.2	514.0	800.0	1,044.3	1,317.7																		
000077E04B-11	3	F	TS	25	118.4	178.8	228.2	234.2	343.9	530.5	943.2	1,195.9	1,344.4	1,809.7																	
000077E080-22	3	F	TS	27	89.2	168.7	158.8	197.6	285.9	492.7	686.7	831.7	1,040.8	1,520.0	1,671.3																
000077E043-16	3	F	TS	27	60.3	202.2	180.3	206.7	251.3	278.4	408.6	467.4	719.1	1,059.5	1,580.5																
000077E0F2-13	3	F	TS	27	88.9	142.2	156.1	206.1	311.3	425.3	595.6	742.6	976.6	1,298.4	1,705.2																
000077E0C4-5	3	F	TS	32	118.6	97.5	104.4	178.5	240.9	367.1	469.3	550.9	679.8	905.9	1,137.8	1,288.0	1,552.1														
000077E0873-26	4	F	TS	25	102.3	175.7	238.2	255.9	411.9	634.4	840.4	1,038.5	1,308.3	1,574.7																	
000077E02A-28	4	F	TS	27	51.6	96.2	130.2	161.2	238.2	320.3	453.4	531.6	800.3	1,261.3	1,671.5																
000077E054-2	4	F	TS	34	109.7	136.7	131.4	205.2	221.8	362.4	422.0	463.7	601.9	839.5	928.5	935.4	1,156.1	1,639.9													
000077E066-23	4	F	TS	25	82.9	169.1	195.3	231.9	332.8	488.6	691.9	862.5	1,055.0	1,302.9																	
000077E0CC-29	4	F	TS	39	100.7	97.7	117.2	150.3	228.3	278.4	402.0	423.8	489.3	550.9	650.5	660.7	810.7	939.4	1,187.6	1,536.6											
000077E0814-1	5	F	SS	50	111.4	109.7	128.4	125.3	153.6	176.4	203.4	185.9	222.9	244.5	248.8	261.6	281.3	274.4	281.5	330.7	220.5	228.1	260.9	218.1	202.2						
000077E0800-27	5	F	TS	41	85.0	188.7	162.7	210.8	263.8	429.3	553.2	711.5	866.2	1,146.9	1,184.1	1,122.3	1,237.1	1,187.8	1,333.7	1,439.4	1,800.0										
000077E0840-19	5	F	TS	48	61.2	115.3	165.5	210.9	291.2	382.2	528.2	593.0	786.5	932.6	903.5	1,049.3	1,145.2	1,070.3	1,079.6	1,225.3	1,266.5	1,321.4	1,376.0	1,539.9							
000077E0889-25	5	F	TS	22	136.5	198.8	296.4	316.2	525.3	752.3	1,039.5	1,185.5	1,668.2																		
000077E0CC-49	5	F	SS	50	82.5	75.7	121.1	72.5	178.6	312.7	330.7	402.7	517.5	588.8	523.4	523.4	440.1	402.7	527.8	538.1	450.0	452.8	423.8	427.2	436.6						

Tumour ulceration grade 1

Tumour ulceration grade 2

Figure 33: Individual group plots





Figure 33(contd): Individual group plots





Figure 33(contd): Individual group plots




Figure 34: Tumour weight

Group	Animal ID	Date	Day	Final Tumour Weight/ g	Endpoint
1	00077E6CC9-7	21/05/2014	27	1.347	Oversized tumour
1	00077E7D9C-3	19/05/2014	25	1.761	Oversized tumour
1	00077E7ED6-17	04/06/2014	41	1.128	Oversized tumour
1	00077E7FDB-14	13/06/2014	50	No tumour present	Study Termination
1	00077E8B08-20	21/05/2014	27	1.197	Oversized tumour
2	00077E6AFE-21	13/06/2014	50	No tumour present	Study Termination
2	00077E6C4D-12	13/06/2014	50	0.091	Study Termination
2	00077E7EF9-18	26/05/2014	32	1.882	Oversized tumour
2	00077E8B5F-30	19/05/2014	25	1.432	Oversized tumour
2	00077E8C0D-4	16/05/2014	22	1.119	Oversized tumour
3	00077E6ABD-11	19/05/2014	25	1.717	Oversized tumour
3	00077E6B0F-22	21/05/2014	27	1.398	Oversized tumour
3	00077E6D43-16	21/05/2014	27	1.098	Oversized tumour
3	00077E7DF2-13	21/05/2014	27	1.45	Oversized tumour
3	00077E8C4C-5	26/05/2014	32	1.185	Oversized tumour
4	00077E6B73-26	19/05/2014	25	1.714	Oversized tumour
4	00077E6D2A-28	21/05/2014	27	1.399	Oversized tumour
4	00077E6D5A-2	28/05/2014	34	1.42	Oversized tumour
4	00077E6D6E-23	19/05/2014	25	1.42	Oversized tumour
4	00077E7DCC-29	02/06/2014	39	1.781	Oversized tumour
5	00077E6B14-1	13/06/2014	50	0.164	Study Termination
5	00077E6BD0-27	04/06/2014	41	1.474	Oversized tumour
5	00077E80AD-19	11/06/2014	48	1.418	Oversized tumour
5	00077E8A89-25	16/05/2014	22	1.143	Oversized tumour
5	00077E8CC4-9	13/06/2014	50	0.297	Study Termination

REFERENCES CITED IN THE DESCRIPTION

This list of references cited by the applicant is for the reader's convenience only. It does not form part of the European patent document. Even though great care has been taken in compiling the references, errors or omissions cannot be excluded and the EPO disclaims all liability in this regard.

Patent documents cited in the description

- WO 2009105044 A [0006] [0215]
- WO 2013150303 A [0007]
- WO 2014013258 A [0007]
- US 5851832 A [0055]
- US 6777233 B [0055]
- US 6468794 B [0055]
- US 5753506 A [0055]
- WO 2005121318 A [0055]
- US 7416888 B [0154] [0159]
- EP 1645626 B1 [0156] [0159]
- US 7514259 B [0165] [0214]
- GB 2013050879 W [0181]
- WO 2012004611 A [0202]
- GB 1314573 A [0438]
- GB 1317887 A [0438]

Non-patent literature cited in the description

- HASSANI ; O'REILLY ; PEARSE ; STROEMER et al. *PLoS One*, 2012, vol. 7 (11) [0003]
- FONSATO et al. *Stem Cells*, September 2012, vol. 30 (9), 1985-98 [0007]
- MILJAN et al. *Stem Cells Dev.*, 2009 [0020]
- KATSUDA et al. *Proteomics*, 2013 [0073]
- KATSUDA et al. *Scientific Reports*, 2013 [0073]
- GENNARO. Remington: The Science and Practice of Pharmacy. 2000 [0169] [0437]
- E W MARTIN. Remington's Pharmaceutical Sciences [0171]
- ISHIWATA et al. *World J Gastroenterol.*, 28 January 2011, vol. 17 (4), 409-418 [0190]
- KAUR et al. *BMC cancer*, 2012, vol. 12, 120 [0190]
- ROFSTAD BR. *J. Cancer*, 1994, vol. 70, 804-812 [0190]
- JACOBS et al. *ASN Neuro.* 2011, 2011, vol. 3 (3) [0190]
- MOKRY et al. *Stem Cells Dev.*, 2004, vol. 13, 658-664 [0191]
- WORKMAN et al. *British Journal of Cancer*, 2010, vol. 102, 1555-1577 [0404]
- AMBROS et al. *RNA*, 2003, vol. 9, 277-279 [0437]
- BANERJEE, S. ; WILLIAMSON, D. ; HABIB, N. ; GORDON, M. ; CHATAWAY, J. *Age and Ageing*, 2011, vol. 40, 7-13 [0437]
- CHUNG et al. *Cell Stem Cell*, 2008, vol. 2, 113-117 [0437]
- DAI, L.J. ; MONIRI, M.R. ; ZENG, Z.R. et al. *Cancer Lett*, 2011, vol. 305 (1), 8-20 [0437]
- DING, D.C. ; SHYU, W.C. ; LINS, Z. *Cell Transplant*, 2011, vol. 20, 5-14 [0437]
- EINSTEIN, O. ; BEN-HUR, T. *Arch Neurol*, 2008, vol. 65, 452-456 [0437]
- HASSANI Z ; O'REILLY J ; PEARSE Y ; STROEMER P ; TANG E ; SINDEN J ; PRICE J ; THURET S. Human neural progenitor cell engraftment increases neurogenesis and microglial recruitment in the brain of rats with stroke. *PLoS One*, 2012, vol. 7 (11), e50444 [0437]
- HODGES et al. *Cell Transplant.*, 2007, vol. 16 (2), 101-15 [0437]
- HORIE, N. ; PEREIRA, N.P. ; NIIZUMA, K. ; SUN, G. et al. *Stem Cells*, 2011, vol. 29, 274-285 [0437]
- HULKOWER, K. I. ; HERBER, R. L. *Pharmaceutics*, 2011, vol. 3, 107-124 [0437]
- KATSUDA ; KOSAKA ; TAKESHITA ; OCHIYA. *Proteomics*, 2013, vol. 00, 1-17 [0437]
- KATSUDA ; TSUCHIYA ; KOSAKA ; YOSHIOKA ; TAKAGAKI ; OKI ; TAKESHITA ; SAKAI ; KURODA ; OCHIYA. *Scientific Reports*, 2013, vol. 3 (1197), 1-11 [0437]
- KORNBLUM. *Stroke*, 2007, vol. 38, 810-816 [0437]
- LAI et al. Proteolytic Potential of the MSC Exosome Proteome: Implications for an Exosome-Mediated Delivery of Therapeutic Proteasome. *International Journal of Proteomics*, 2012, 14 [0437]
- LITTLEWOOD, T. D. ; HANCOCK, D. C. ; DANIELIAN, P. S. et al. *Nucleic Acid Research*, 1995, vol. 23, 1686-1690 [0437]
- MILJAN, E.A. ; SINDEN, J.D. *Current Opinion in Molecular Therapeutics*, 2009, vol. 4, 394-403 [0437]
- MILJAN EA ; HINES SJ ; PANDE P ; CORTELING RL ; HICKS C ; ZBARSKY V ; UMACHANDRAN, M ; SOWINSKI P ; RICHARDSON S ; TANG E. Implantation of c-mycER TAM immortalized human mesencephalic-derived clonal cell lines ameliorates behavior dysfunction in a rat model of Parkinson's disease. *Stem Cells Dev.*, March 2009, vol. 18 (2), 307-19 [0437]

- MITCHELL et al. *Journal of Immunological Methods*, 2008, vol. 335, 98-105 [0437]
- POLLOCK et al. *Exp Neurol.*, May 2006, vol. 199 (1), 143-55 [0437]
- MARK F PITTENGER ; ALASTAIR M MACKAY ; STEPHEN C BECK ; RAMA K JAISWAL et al. *Science*, 02 April 1999, vol. 284, 5411 [0437]
- SHAH, K. *Adv Drug Deliv Rev*, 2012, vol. 64 (8), 739-748 [0437]
- SMITH, E. J. ; STROEMER, R.P. ; GORENKOVA, N. ; NAKAJIMA, M. et al. *Stem Cells*, 2012, vol. 30, 785-796 [0437]
- STEVENATO, L. ; CORTELING, R. ; STROEMER, P. ; HOPE, A. et al. *BMC Neuroscience*, 2009, vol. 10, 86 [0437]
- STROEMER, P. ; PATEL, S. ; HOPE, A. ; OLIVEIRA, C. ; POLLOCK, K. ; SINDEN, J. *Neurorehabil Neural Repair*, 2009, vol. 23, 895-909 [0437]
- THÉRY, C. ; OSTROWSKI, M. ; SEGURA, E. et al. *Nature Reviews Immunology*, 2009, vol. 9, 581-593 [0437]
- THEIR et al. Direct Conversion of Fibroblasts into Stably Expandable Neural Stem Cells. *Cell Stem Cell*, 20 March 2012 [0437]
- TIMMERS, L. ; LIM, S. K ; ARSLAN, F. ; ARMSTRONG, J. S. et al. *Stem Cell Res*, 2007, vol. 1, 129-137 [0437]
- YUAN, S.J. ; MARTIN, J ; ELIA, J. ; FLIPPIN, J. et al. *PLoS ONE*, 2011, vol. 6, e17540 [0437]